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Tokuyama Dental

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016 JAPAN
Phone: +81 3-3835-2261 Fax: +81 3-3835-2265

27th of May, 2020
«Approved by»

Makoto Sakamoto

President

Tokuyama Dental Corporation

Instruction for Use

Dental composite material OMNICHROMA

Versions:

1. Dental composite material OMNICHROMA consisting of:
 - 1.1. Dental composite material OMNICHROMA in syringes 4g / 2.2 ml;
 - 1.2. Instruction leaflet – 1pc.
2. Dental composite material OMNICHROMA in unidoses consisting of:
 - 2.1. Dental composite material OMNICHROMA in unidoses (PLT) 0.2 g (20 pcs.)
 - 2.2. Instruction leaflet – 1pc.

Revision 2

1. Name of the medical device.

Dental composite material OMNICHROMA

Versions:

1. Dental composite material OMNICHROMA consisting of:
 - 1.1. Dental composite material OMNICHROMA in syringes 4g / 2.2 ml;
 - 1.2. Instruction leaflet – 1pc.
2. Dental composite material OMNICHROMA in unidoses consisting of:
 - 2.1. Dental composite material OMNICHROMA in unidoses (PLT) 0.2 g (20 pcs.)
 - 2.2. Instruction leaflet – 1pc.

Hereinafter referred as – the material, device, material OMNICHROMA, OMNICHROMA

2. Details of the designer, manufacturer and the authorized representative

Designer

Tokuyama Dental Corporation
38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan
Tel: +81-3-3835-7201
URL: <http://www.tokuyama-dental.com>

Manufacturer

Tokuyama Dental Corporation
38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan
Tel: +81-3-3835-7201
URL: <http://www.tokuyama-dental.com>

Manufacturing site

Tokuyama Dental Corporation Kashima Factory
26, Sunayama, Kamisu-shi, Ibaraki 314-0255, Japan

Authorized representative of the manufacturer in the Russian Federation

«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter “I”

Operating address: 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"

Phone +7 (812) 635-88-90; Email: info@protecodent.ru

3. Intended use of the medical device with indication of the potential consumer, field of application.

Intended use: This product is suitable for conventional restoration of tooth cavities

Potential consumers of the medical device: dentists.

Field of application: dentistry.

4. Indications, contraindications and possible side effects associated with the intended use of the medical device.

Indications:

For use as a tooth shade resin material in dental procedures, such as:

- Direct anterior and posterior restorations
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

Contraindications:

OMNICHROMA contains methacrylic monomers and UV absorbers. Do not use OMNICHROMA for patients allergic to or hypersensitive to methacrylic and related monomers, UV absorbers or any of the other ingredients.

Possible side effects: allergic reactions to or hypersensitive to methacrylic and related monomers, UV absorbers or any of the other ingredients.

5. Description and operating principal of the medical device.

Dental composite material OMNICHROMA is a light-cured, radiopaque composite resin. OMNICHROMA is a single shade material.

OMNICHROMA contains 79% by weight (68% by volume) of spherical silica-zirconia filler (mean particle size: 0.3 μm , particle size range: 0.2 to 0.4 μm) and composite filler. Monomer matrix of OMNICHROMA contains 1,6-bis(methacryl-ethyloxycarbonylamino) trimethyl hexane (UDMA), Triethylene glycol dimethacrylate (TEGDMA), Mequinol, Dibutyl hydroxyl toluene and UV absorber. It is recommended to use in anterior and posterior restorations and is indicated for all carious classes.

OMNICHROMA is available in either syringes 4g/2,2 ml or in unidoses "Pre Loaded Tip" (PLT) 0,2 g.

6. Precautionary measures

- 1) DO NOT use OMNICHROMA for any purpose other than those listed in these instructions. Use OMNICHROMA only as directed herein.
- 2) The material is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- 3) DO NOT use OMNICHROMA if the safety seals are broken or appear to have been tampered with.
- 4) If OMNICHROMA causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling OMNICHROMA to avoid the possibility of allergic reactions from methacrylic monomers. Notes: Certain substance/materials may penetrate through examination

gloves. If OMNICHROMA comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.

- 6) Avoid contact of OMNICHROMA with eyes, mucosal membrane, skin and clothing.
- If OMNICHROMA comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
- If OMNICHROMA comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
- If OMNICHROMA comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
- Instruct the patient to rinse his mouth immediately after treatment.
- 7) OMNICHROMA should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 8) To avoid the unintentional ingestion of OMNICHROMA, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments and brushes with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 11) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.
- 12) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
- 13) DO NOT mix OMNICHROMA with other brands of resin composite to avoid losing its color adjustment.

7. Order of application

7.1. Special notes for the use of unidoses (PLT- pre-loaded tip)

- 1) Unidoses (PLTs) are designed for single patient use only. Do not re-cap and/or re-use the unidose tip once material has been dispensed for that patient.
- 2) Dispensers are not provided for OMNICHROMA. Use a dispenser that is permitted for use in accordance with the standard procedure and fits the unidose of OMNICHROMA. For proper use and proper disinfection, see manufacturer's instructions.
- 3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

7.2. Dental composite material OMNICHROMA in syringes 4 g / 2.2 ml and unidoses (PLT) 0.2 g.

Determine if this product is suitable for application to the case before use. The desirable esthetic outcome may not be obtained depending on the case (e.g. a case in which no tooth structure remains on the lingual surface, a case in which the underlying tooth is extremely discolored, a case in which the color of the underlying tooth is to be changed to another color, etc.).

1. Cleaning

Instruction for Use of the medical device
Dental composite material OMNICHROMA

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Isolation

A rubber dam is the preferred method of isolation.

3. Cavity Preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II), bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- Scalloped bevel could be desirable depending on the esthetics and retention.
- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

4. Pulp Protection

Glass ionomer or a calcium hydroxide lining should be applied if the cavity is in close proximity to the pulp. **DO NOT USE EUGENOL-BASED MATERIALS** to protect the pulp as these materials will inhibit curing OMNICHROMA. OMNICHROMA is not indicated for direct pulp capping.

5. Bonding System

Apply a bonding system according to its manufacturer's instructions.

If you want to use a self-cured or dual-cured bonding system, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

6. Dispensing

Unidoses (PLT)

Please read **SPECIAL NOTES FOR THE USE OF** unidoses (PLT) before use.

- Load the PLT into a dispenser that fits the PLT
- Remove the PLT cap.
- Extrude the paste directly into the cavity or indirectly on the mixing pad.

Syringe

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.

- **DO NOT** apply unnecessary force to the syringe immediately after removing from the refrigerator.

7. Filling and Contouring

Fill the cavity incrementally. Increments should not exceed the indicated curing depth (please refer to Table 1).

- Use OMNICHROMA BLOCKER under OMNICHROMA for the direct restoration of extensive Class III or IV cavity. OMNICHROMA can appear semi-translucent when OMNICHROMA is not adequately surrounded by a cavity wall. If applied in direct restorations of extensive Class III or IV cavity, the cured OMNICHROMA may appear dark.

- In the case of repair of porcelain / composite, the application of an opaquer to the fractured surfaces as an initial layer is needed to hide a metal portion if it exists. Then, place OMNICHROMA into the fractured surfaces.

- In the case of masking tooth discoloration, apply OMNICHROMA BLOCKER or a dental color masking material such as an opaquer on the discolored surface to mask the discoloration, then place OMNICHROMA.

8. Curing

Light-cure each increment for at least the indicated time (please see Table 1), keeping the curing light tip within a distance of 2 mm from the increment.

9. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

7.3. Indications for filling and curing

The material is designed to be used with either a halogen or LED curing-light with a wavelength of 400-500 nm. Be sure to light-cure the material extra-orally and check the time needed for complete hardening of the material with your light-curing unit before performing the bonding procedure. The below table illustrates the relationship between curing time and increment depth.

Table 1. Relationship between curing time and increment depth of OMNICHROMA.

Light type	Intensity (mW/cm2)	Curing time (seconds)	Increment depth (mm) ¹
Halogen	400	20	2.7
		40	3.3
	600	20	3.1
	800	10	2.7
		20	3.3
LED	400	20	2.5

		40	3.6
	600	20	2.9
	800	10	2.6
		20	3.1
	1200	10	3.0
	2300	6	2.9
¹ Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO 4049.			

8. Performance specifications.

Table 2. Dental composite material OMNICHROMA.

Test item	Specification ¹⁾
Sensitivity to ambient light	The material shall remain physically homogeneous.
Depth of cure	Minimum 1.5mm
Flexural strength	Minimum 80 MPa
Water sorption	Maximum 40 µg/mm ³
Solubility	Maximum 7.5 µg/mm ³
Color stability	No more than a slight change in color shall be observed.
Radio-opacity	The radio-opacity shall be equal to or greater than of the same thickness of aluminum.

¹⁾ Certificates of Analysis (see Annex 3)


The above tests are performed by referring to ISO 4049. The Film thickness of ISO 4049 is not applied as a test item for OMNICHROMA because this product is not a luting material. The Working time and Setting time of ISO 4049 are not applied as test items for OMNICHROMA because this product does not belong to the Class 1 or Class 3 material defined in ISO 4049. The Shade of ISO 4049 is not applied as a test item for OMNICHROMA because this product is intended to be supplied not in various shades.

Dental composite material OMNICHROMA is functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

9. The standards of compliance of the medical device

ISO 13485:2016, ISO 14971:2012, ISO 15223-1:2016, ISO 7405:2008, EN ISO 10993-1-2009, EN ISO10993-3-2014, EN ISO 10993-5-2009, ISO 10993-10:2010, EN ISO 10993-12-2012, ISO 4049:2009, EN ISO 10993-11-2009, EN ISO 10993-18-2009, EN 1041-2008, EN 1641-2009

10. Interpretation of symbols used in labeling of consumer and transport packaging:

Symbol	Description	Symbol	Description
	Manufacturer		Lot number
	Catalogue number		"Use by date..."
	Warning symbol to show the possibility for the occurrence of skin allergic reaction		«Caution! Refer to the instructions for use»
	Temperature range		Authorized Representative in the EU
	Compliance with European standards		mark of conformity GOST R

11. Requirements for maintenance and repair of the medical device.

Medical devices do not required maintenance and repair.

12. Sterilization, cleaning and disinfection.

The medical device is supplied non-sterile and is not intended to sterilize the product before using it.

13. Conditions of transportation, storage and operation.

Transportation of products packaged in shipping containers must be carried out by all types of closed vehicles in compliance with the precautions specified on the shipping container and in accordance with the rules for the carriage of goods in force on this type of transport.

Transportation terms and conditions:

- temperature - from 0 to + 25°C;
- relative humidity of air - 10- 100%
- Atmosphere pressure: 800 ~ 1,060 hPa

Storage conditions:

- temperature - from 0 to + 25°C;
- relative humidity of air - 10-70 %
- Atmosphere pressure: 900 ~ 1,060 hPa

Operating Conditions

- Operation temperature from +32 ° C to +42 ° C

Room conditions during operation of the medical devices:

- Temperature +18 °C-+30 °C

14. Shelf life

Shelf life: 3 years from a date of production.

15. The procedure for the disposal and destruction of medical devices

Medical device Dental composite material OMNICHROMA, depending on the degree of their epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe or unidose (PLT) should be removed and light-cured before disposing.

The products must be collected and transferred to specialized licensed organizations involved in the disposal of medical waste. Transportation of garbage from institutions to the place of disposal is carried out by specialized companies.

Medical devices must be disposed in accordance with medical procedures in medical institutions and local, state or federal regulations. The choice of methods for safe disposal and destruction of products depends on the capacity and profile of the medical organization, the availability of facilities for decontamination / neutralization of waste, method of disposal / destruction of waste, adopted in the administrative territory (incineration, landfill, disposal).

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This public document

2. has been signed by HATANO Ryuji

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of HATANO Ryuji ,Notary

Certified

5. at Tokyo

6. JUN. 03. 2020

7. by the Ministry of Foreign Affairs

8. 20- N2016338

9. Seal/stamp:

10. Signature

Tanaka Toshie

TANAKA Toshie

For the Minister for Foreign Affairs

Registered No. 260 of 2020

Certificate of Acknowledgment of Notary

On this 3rd day of June , 2020 , before me , HATANO
Ryuji , a notary in and for Tokyo Legal Affairs Bureau ,
personally appeared Tomoko MIFUNE, with satisfactory
evidence of her identification and of her being the agent of
Makoto SAKAMOTO, President of Tokuyama Dental
Corporation, said Makoto SAKAMOTO , had acknowledged
to have affixed his signature to the attached document.
Witness , I set my hand and seal .

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Hereinafter referred as – the material, device, material OMNICHROMA BLOCKER, OMNICHROMA BLOCKER.

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«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I"

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Indications:

For use as a tooth shade resin in dental procedures.

OMNICHROMA BLOCKER contains methacrylic monomers. Do not use OMNICHROMA BLOCKER for patients allergic to or hypersensitive to methacrylic and related monomers or any other ingredients.

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Dental composite material OMNICHROMA BLOCKER is a light-cured, radiopaque, composite resin. OMNICHROMA BLOCKER contains 82% by weight (71% by volume) of silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. All inorganic filler contained in OMNICHROMA BLOCKER is a spherical filler (mean particle size: 0.2 μm , particle size range: 0.1 to 0.3 μm). The monomer matrix contains Bis-GMA and Triethylene glycol dimethacrylate.

OMNICHROMA BLOCKER is designed for use as a thin layer at the lingual cavity wall of extensive class III and IV restorations. OMNICHROMA BLOCKER can mask slight discoloration or reconstruct a highly opaque tooth.

OMNICHROMA BLOCKER is provided either in syringe 4g/2.1 ml or in unidoses "Pre Loaded Tip" (PLT) 0,2 g.

6. Precautionary measures

- 1) DO NOT use OMNICHROMA BLOCKER for any purpose other than those listed in these instructions. Use OMNICHROMA BLOCKER only as directed herein.
- 2) The material is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- 3) DO NOT use OMNICHROMA BLOCKER if the safety seals are broken or appear to have been tampered with.
- 4) If OMNICHROMA BLOCKER causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling OMNICHROMA BLOCKER to avoid the possibility of allergic reactions from methacrylic monomers. Notes: Certain substance/materials may penetrate through examination gloves. If OMNICHROMA BLOCKER comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.

- 6) Avoid contact of OMNICHROMA BLOCKER with eyes, mucosal membrane, skin and clothing.
 - If OMNICHROMA BLOCKER comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If OMNICHROMA BLOCKER comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If OMNICHROMA BLOCKER comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) OMNICHROMA BLOCKER should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 8) To avoid the unintentional ingestion of OMNICHROMA BLOCKER, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments and brushes with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 11) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application

7.1. Special notes for the use of unidoses (PLT)

PLT is a pre-loaded tip.

- 1) Unidoses (PLTs) are designed for single patient use only. Do not re-cap and/or re-use the unidose tip once material has been dispensed for that patient.
- 2) Dispensers are not provided for OMNICHROMA BLOCKER. Use a dispenser that is permitted for use in accordance with the standard procedure and fits the unidose of OMNICHROMA BLOCKER. For proper use and proper disinfection, see manufacturer's instructions.
- 3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

7.2. Dental composite material OMNICHROMA BLOCKER in syringes 4g / 2.1 ml and in unidoses (PLT) 0.2 g (20 pcs.).

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2. Shade Selection

Select OMNICHROMA BLOCKER as the initial layer under OMNICHROMA, with direct restorations of extensive Class III or IV cavity using OMNICHROMA.

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- Scalloped bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

5. Pulp Protection

Glass ionomer or a calcium hydroxide lining should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of bonding agent and/or OMNICHROMA BLOCKER. OMNICHROMA BLOCKER is not indicated for direct pulp capping.

6. Bonding System

Apply a bonding system according to manufacturer's instructions.

- If you want to use a self-cured or dual-cured bonding system, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

7. Dispensing

Unidoses (PLT)

Please read SPECIAL NOTES FOR THE USE OF PLT before use.

- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.
- Extrude the paste directly into cavity or indirectly on the mixing pad.

Syringe

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately after dispensing.

- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling

Fill the cavity with OMNICHROMA BLOCKER as an initial layer. Increments should not exceed the indicated depth (please refer to Table 1).

- DO NOT mix OMNICHROMA BLOCKER with other brands of resin composite to avoid incomplete cure or voids.

9. Curing

Light-cure each increment for at least the indicated time (please see Table 1), keeping the curing-light tip within a distance of 2 mm from the increment.

- Follow the filling and curing instructions of OMNICHROMA to be layered over the cured OMNICHROMA BLOCKER.

10. Finishing

Shape and polish the restoration according to instructions for OMNICHROMA used as the surface layer of the restoration.

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Table 1. Relationship between curing time and increment depth of OMNICHROMA BLOCKER.

	Intensity	Curing time	Increment depth
--	-----------	-------------	-----------------

Light type	(mW/cm ²)	(seconds)	(mm) ¹
Halogen	400	20	1.8
	600	20	1.9
	800	10	1.7
		20	2.0
LED	400	20	1.6
	600	20	1.7
	800	10	1.7
		20	1.9
	1200	10	1.8
	2300	6	1.9

¹ Increment depth was determined on the basis of test results performed in accordance with section 7.10 “depth of cure” of ISO 4049.

8. Performance specifications.

Table 2. Dental composite material OMNICHROMA BLOCKER

Test item	Specification ⁴⁾
Sensitivity to ambient light ²⁾	After light exposure for 60 ± 5 sec, the material of samples remains physically homogeneous
Depth of cure ²⁾	Minimum 1.0mm (Opaque shades)
Flexural strength ²⁾	Minimum: 80MPa
Water sorption ²⁾	Maximum: 40µg/mm ³
Solubility ²⁾	Maximum: 7.5µg/mm ³
Shade ²⁾	The shade of the set material which is immersed into water for 7 days shall match closely that of the manufacturer's shade guide.
Appearance ³⁾	The set material shall be homologous without bubbles and contaminations

Color stability ¹⁾	No more than a slight change in color shall be observed after light irradiation.
Radio-opacity ²⁾	The radio-opacity of the set material shall be equal to or greater than that of the same thickness of aluminium

1) This test was performed according to the standard ISO 4049.

2) This test was performed according to the standard ISO 4049.

3) These tests were performed according to the internal standards of the manufacturer

4) Certificates of Analysis (see Annex 3)

Dental composite material OMNICHROMA BLOCKER is functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

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Storage conditions:

- temperature - from 0 to + 25°C;
- relative humidity of air - 10-70 %

Instruction for Use of the medical device
Dental composite material OMNICHROMA BLOCKER

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
Medical device Dental composite material OMNICHROMA BLOCKER, depending on the degree of it's epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe or unidose (PLT) should be removed and light-cured before disposing.

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4. bears the seal/stamp of HATANO Ryuji , Notary

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Notary, HATANO Ryuji

Ueno Notary's Office of Tokyo Legal
Affairs Bureau , 1-7-2 Higashiueno ,
Taito-ku , Tokyo, Japan

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Tokuyama Dental Corporation

Supplement No. 1

To the Extract from the technical documentation

Dental composite material OMNICHROMA

1. Clause 1. «Name of medical device» shall be amended to read as follows:

Dental composite material OMNICHROMA

Versions:

1. Dental composite material OMNICHROMA consisting of:
 - 1.1. Dental composite material OMNICHROMA in syringes 4g / 2.2 ml;
 - 1.2. Instruction leaflet – 1pc.
2. Dental composite material OMNICHROMA in unidoses consisting of:
 - 2.1. Dental composite material OMNICHROMA in unidoses (PLT) 0.2 g (20 pcs.)
 - 2.2. Instruction leaflet – 1pc.
3. Dental composite material OMNICHROMA BLOCKER consisting of:
 - 3.1. Dental composite material OMNICHROMA BLOCKER in syringes 4g / 2.1 ml;
 - 3.2. Instruction leaflet – 1 pc.
4. Dental composite material OMNICHROMA BLOCKER in unidoses consisting of:
 - 4.1. Dental composite material OMNICHROMA BLOCKER in unidoses (PLTs) 0.2 g (20 pcs.).
 - 4.2. Instruction leaflet – 1 pc.

2. Clause 14.3 «Draft labeling» shall be supplemented by drafts of labels for primary, secondary and transport packages and shall be amended to read as follows:

Primary packaging labeling (syringes)

- The shortened name of a medical device (allowed identifying the product if a sticker has limited size).
- Name and country of the manufacturer.
- Name of the manufacturer representative in the EU and the USA (if applicable)
- Volume and / or mass of contents of a syringe
- Lot number (LOT) of a primary packaging, syringe with material (differs from a lot number (LOT) of a secondary (consumer) package.
- Catalogue number (REF)
- Storage conditions (Temperature range).
- CE mark (if applicable)
- Expiration Date («Use before» Symbol)
- Symbol - «Caution! Refer to the instructions for use »

Label draft

OMNICHROMA		4 g (2.2 mL)
Tokuyama Dental Corporation, Japan		
LOT	 ГОД, МЕСЯЦ	REF 
(0 до +25°C)		

OMNICHROMA BLOCKER		4 g (2.1 mL)
Tokuyama Dental Corporation, Japan		
LOT	 год, месяц	REF 
(0 до +25°C)		

Primary packaging labeling (unidoses)

- The shortened name of a medical device (allowed to identify the product if a sticker has limited size).
- Lot number (LOT) of a primary packaging, syringe with material (differs from a lot number (LOT) of a secondary (consumer) package)
- Expiration Date («Use before» Symbol)

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
Marking of secondary (consumer) packaging (syringes):

- Name of a Medical Device.
- Name and address of the manufacturer.
- Name and address of a representative of the EU or the USA.(if applicable).
- Name and address of the authorized representative of the manufacturer in the Russian Federation.
- Storage conditions (symbol Temperature range).
- Material classification according to ISO 4049, radiopacity
- Mass, g and / or volume, ml of the contents of the syringe.
- The expiry date of the material (symbol “Use to”).
- The presence of hazardous substances: symbols with the text of caution. (if applicable)
- Lot number (LOT) of a primary packaging, syringe with material (differs from a lot number (LOT) of a secondary (consumer) package).
- CE mark (if applicable)
- Symbol - “Caution! Refer to the instructions for use »

- Number and date of a registration certificate.
- PCT mark (symbol).

Label draft


Marking of secondary (consumer) packaging (unidoses):

- Name of a Medical Device.
- Name and address of the manufacturer
- Name and address of a representative of the EU or the USA.(if applicable)
- Name and address of the authorized representative of the manufacturer in the Russian Federation.
- Lot number (LOT) of a primary packaging, syringe with material (differs from a lot number (LOT) of a secondary (consumer) package).
- The expiry date of the material (symbol "Use to").
- CE mark (if applicable)
- Mass, g of the contents of a unidoses (PLT), quantity.
- Storage conditions (symbol Temperature range).
- Material classification according to ISO 4049, radiopacity
- The presence of hazardous substances: symbols with the text of caution. (if applicable)
- Symbol - "Caution! Refer to the instructions for use »
- Number and date of a registration certificate.
- PCT mark (symbol).

Marking of transport packaging:

- Name of a Medical Device (Not labeled if various medical devices are contained in one transport package.)
- Manufacturer name and address (city and country)

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APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by HATANO Ryuji
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of HATANO Ryuji ,Notary
5. at Tokyo
6. JUN. 03. 2020
7. by the Ministry of Foreign Affairs
8. 20- N0016335
9. Seal/stamp:
10. Signature

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