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«APPROVED»

«B. Well Swiss AG»

Bernadette Maeder-Holenstein

B. Maeder-Holenstein

INSTRUCTION MANUAL

Blood pressure measuring devices models PRO-60, MED-61,

MED-62, MED-63 with accessories

[Handwritten mark]

1. INTRODUCTION

Thank you for purchasing a B. Well blood pressure monitor! The B. Well mechanical tonometer is reliable, durable and provides high measurement accuracy. Blood pressure measuring devices models PRO-60/MED-61/MED-62/MED-63 is designed for indirect measurement of systolic and diastolic blood pressure by the Korotkov method.

Area of application: Medical and preventive health care institutions and individual use for self-monitoring.

Important advantages of PRO-60 / MED-61 / MED-62 / MED-63:

- Durable and comfortable cuff.
- High-precision manometer in metal casing.
- Valve with smooth air outlet.

2. CLASSIFICATION OF ARTERIAL PRESSURE MEANS

Table for the classification of blood pressure values (unit: mmHg) in accordance with the European Society of Hypertension (ESH)

Range	Systolic pressure blood	Diastolic pressure blood	Actions
3 degree: severe form of essential hypertension	180 or higher	110 или выше	Urgently consult a doctor!
2 degree: the middle form of hypertensive disease	160-179	100-109	Get medical attention immediately
1 degree: mild form of hypertension	140-159	90-99	Contact your doctor
The upper limit of the norm	130-139	85-89	Contact your doctor
Normally	Below 130	Below 85	Self-monitoring
Optimal	Below 120	Below 80	Self-monitoring

① NOTE: Show the measured values to your doctor. Never use the results of your measurements to independently change the doses of drugs prescribed by your doctor.

3. COMPOSITE PARTS OF THE APPLIANCE

The stethoscope for tonometers is delivered depending on the equipment. The cuff and, accordingly, the pneumatic chamber, can be executed in different sizes depending on a complete set.

In the MED-63 model, the head of the stethoscope is built into the cuff.

Choice of cuff size.

Choosing the right cuff size is an important condition that affects the accuracy of the measurements. The cuff size should correspond to the circumference of your shoulder (hip), measured with a tight fit in the middle of the shoulder (hip).

The table shows the correspondence of the sizes of the cuffs and the girth of the shoulder (thigh).

1	2	3	4	5	6	7	8	9
Cuff name	N	I	C	S	M	M-L	L	T
Girth of the shoulder (hip), cm	9-20	11-22	14-25	18-30	22-32	22-42	32-52	40-66

4. PRECAUTIONS

This appliance is intended for use solely for the purposes specified in this manual.


The manufacturer is not liable for damages caused by improper use.

Never inflate the cuff to values exceeding 300 mm hg. art.

Take care that children do not use the appliance unattended; the device contains small parts.

5. ASSEMBLY INSTRUCTIONS

Connect the two tubes coming from the cuff, respectively the supercharger and the manometer, as shown in the figure.

To assemble the stethoscope, attach the stethoscope head to the single end of the Y-tube (in the MED-63 model, the head of the stethoscope is built into the cuff). To the double end of the Y-tube, attach the metal ear tubes of the stethoscope.


① NOTE: the position of the pressure gauge needle not in the center of the zero area indicated by the symbol "□" is not a fault of the device. The gauge needle must be within the zero range that is set by the manufacturer when the instrument is adjusted in the required position to ensure the required accuracy.

6. MEASUREMENT IMPLEMENTATION

6.1 Before starting measurements

Immediately before measuring blood pressure, avoid eating, smoking and any form of physical stress. All these factors affect the measurement results. Try to find time and relax, sitting in a chair in a relaxed atmosphere for a few minutes before measuring. It is advisable to release the hand that you will use to apply the cuff from the clothing. Do not roll up the sleeve, as it can squeeze your hand, which will lead to inaccuracy in the measurement. Always make measurements on the same hand (usually the left one). Perform measurements regularly at the same time of day, as blood pressure changes throughout the day.


① Warning: use only clinically approved original cuff. A loose cuff can also affect the accuracy of the measurement. Repeated blood pressure measurement should be performed only after a


five-minute break, as blood in the hand stagnates during the measurement, which can lead to incorrect measurement results. Pause depends on age and in elderly people can reach 10-15 minutes.

6.2 Overlay of the cuff

a) Place the cuff on the arm in such a way that its edge is 2-3 cm above the elbow fold.

б) Check that the air tubes go out towards your palm, and the "ARTERY MARK" sign is above the ulnar fossa where the artery passes. Model MED-63: the built-in head of the stethoscope should be 3-5 cm above the elbow bend. Fix the cuff tightly, but not tight. When fixing the cuff on the hand, note whether the mark (white edge or white arrow) falls in the "OK" range, this confirms that the cuff size corresponds to the size of the hand.

в) Put your hand on the table. Make sure that the tubes do not twist.

г) Sit quietly for a few minutes before measuring.

⚠ IMPORTANT! Measurements can be carried out not only on the left, but also on the right hand. In any case, all measurements must be carried out on one hand.

6.3. MEASUREMENT PROCESS

6.3.1 Setting the stethoscope head

For PRO-60 / MED-61 / MED-62 models: place the stethoscope head under the cuff so that it is either below it or 1-2 cm below the cuff. Do not place the stethoscope head on or inside the cuff. Make sure that the head of the stethoscope is in contact with the skin.

For MED-63: Since the head of the stethoscope is built into the cuff, check that it is located 3-5 cm above the elbow bend.

Before using the stethoscope, make sure there are no cracks in the membrane, headphones and tube. Incorrect insertion or damage to the stethoscope will cause distortion of the tone or poor transmission, which will lead to inaccurate measurements.

During the measurement, the body should be pleasantly relaxed. Do not strain your hand during measurement. Make sure that the entry point of the air tube into the cuff is at the heart level.

6.3.2 Pumping the cuff

Close the valve located on the supercharger (pear) by turning the screw clockwise. Do not over tighten. To pump air into the cuff, squeeze the supercharger in your hand evenly, and listen to the pulse at the same time and observe the gauge reading. After you have stopped hearing the pulse, continue to pump air into the cuff, increasing the pressure in it by another 30-40 mm hg. art.

6.3.3 Measurement of systolic blood pressure.

Slowly open the air valve by turning the screw counterclockwise so that the pressure in the cuff falls at a speed of 2-3 mm. gt; Art. (1-2 divisions per manometer) per second. This speed is optimal for obtaining accurate results. When the cuff starts to bleed, you should carefully listen to the tones through the stethoscope. As soon as you hear weak pulse beats in your headphones, remember the manometer readings. This is your systolic (upper) blood pressure value.

6.3.4 Measurement of diastolic blood pressure.

Continue venting the air from the cuff with the same speed (2-3 mm of mercury per second). In the headphones you will hear a pulse, the strength of which changes during the measurement. By the end of the measurement, the strokes will become soft. Watch the pressure gauge needle, at the moment when the sound ceases to be heard, it will show the value of diastolic (lower) blood pressure. After the diastolic pressure value is determined by you, completely release the air from the cuff, remove it from your hand and remove the stethoscope headphones from the ears.

6.3.5 Recording measurements

Repeat the measurement twice. Do not forget to record the values of your measurements, as well as the date and time. The most suitable time for measuring blood pressure is the morning (immediately after sleep), as well as the evening (before dinner). When you visit a doctor, take your notes, this will help him keep your blood pressure under control. Remember that only a doctor is qualified to interpret your pressure.

7. POSSIBLE MALFUNCTIONS AND THEIR METHODS

If there are any problems during the measurement, check the table points and take appropriate measures if necessary.

Poor tone transmission, distortion or extraneous noise

1. Check if the headphones are clogged and if they are not cracked. If not, make sure that they fit tightly and are not worn.
2. Check if the tube has a crack and is not twisted.
3. Check for cracks on the lid and membrane in the working area of the stethoscope.

	4. Make sure that the working part of the stethoscope is in proper contact with the skin and is located above the brachial artery during the measurement. To avoid inaccurate measurements, clean or replace faulty parts.
When inflating the cuff by the supercharger, the pressure does not increase	1. Make sure that the valve is closed. 2. Make sure that the cuff is correctly connected to the supercharger and pressure gauge. 3. Check to see if the cuff, tube and blower have leaks, replace the faulty parts if a malfunction is found.
The air release rate can not be set to 2-3 mm. gt; post by adjusting the air outlet	1. Disconnect the valve from the pear in order to check if there are any obstructions to the air inside the valve. Remove the obstruction and try again. If the valve does not work properly, replace it to avoid inaccurate measurement results.
At rest the pointer is not at 0 ± 3 mm. pillars	1. Make sure that the valve is fully open when checking the zero setting. 2. If the deviation from the zero value exceeds 3 mm. gt; post, contact the sales organization to recalibrate the pressure gauge.

8. CARE OF THE APPLIANCE. CLEANING AND DISINFECTION.

BACKGROUND INFORMATION

With proper care and maintenance, this unit will last you a long time. Follow the general rules described below:

- Do not drop the appliance.
- Do not expose the device to extreme high temperatures, low temperatures, humidity or direct sunlight. Do not touch the fabric from which the cuff is made with sharp tools, as they can damage it.
- Store the cuff, completely releasing the air from it.
- Never disassemble the manometer.
- Store the entire device in a storage bag to keep all parts clean.
- Clean the manometer and the blower with a soft cloth.
- For disinfection, wipe the component parts and accessories with a soft cloth soaked in ethyl alcohol (75-90%) or 3% hydrogen peroxide solution, then dry the cuff in the air. The cuff should only be disinfected with the inner surface of the cuff to prevent damage to the labeling.

Verification information

Verification of instruments for measuring blood pressure is carried out according to document P 50.2.032-2004 "Recommendations on metrology. GSI. Measurers of blood pressure noninvasive. The verification procedure ». Interval between verification: 2 years. Confirmation of the initial verification - an electronic version of the verification certificate, you can find on the site www.bwell-swiss.ru or www.alpha-medica.ru by name, model and serial number of the device. The original of the certificate is in AO Alfa-Medica.

9. MAIN TECHNICAL CHARACTERISTICS

Product Name: Blood pressure measuring devices models PRO-60, MED-61, MED-62, MED-63 with accessories

Pressure measuring range: 20 to 300 mm Hg

Minimum measuring step: 2 mm Hg

Pressure measurement accuracy: ± 3 mmHg

Leakage of air $\leq \pm 3.75$ mm Hg. Art.

Storage temperature: -5 to $+40$ °C

Storage humidity: relative humidity of 15 to 85%

Operating temperature: 0 to 40 °C

Operating humidity: relative humidity from 15 to 85%

Models PRO-60 / MED-61 / MED-62 can be supplied with or without stethoscopes. A cuff without a ring of one of the listed sizes can be delivered to the devices: N, I, C, S, M, M-L, L, T or a cuff with a ring of one of the listed sizes: S, M, M-L, L.

Weight, model PRO-60 / MED-61 / MED-62 / MED-63, not more than: 557.7 g (depending on the type and size of the attached cuff);

The manufacturer reserves the right to change the technical specifications of the device.

10. CONFORMING TO THE STANDARDS. CONTRAINDICATIONS. UTILIZATION.

The high quality of the device is documented:

- Registration Certificate No. RZN YYYY / XXXX from DD.MM.YYYY
- Declaration of conformity.
- Certificate of approval of the type of measuring instruments of the Federal Agency for Technical Regulation and Metrology.

Devices for measuring blood pressure for blood pressure measurement of PRO-60 / MED-61 / MED-62 / MED-63 with accessories meet the following international standards:

- EN 1060-1: 1996 Sphygmomanometers (blood pressure meters) are non-invasive. Part 1. General requirements.
- EN ISO 9001: 2008 "Quality management systems - Requirements"
- EN ISO 13485: 2012 + AC: 2012 "Medical products - Quality management systems - System requirements for regulatory purposes"

Contraindications: no contraindications.

Disposal: The appliance can be disposed of with household waste.

Date of manufacture:

Date of production - the first four digits of the serial number of the device - is indicated on the manometer. The first and second digits are the month number, the third and fourth are the year of production.

11. WARRANTY LIABILITIES

Warranty period for the manometer is 5 years and free service is 10 years. Warranty period for the supercharger, air release valve and stethoscope - 2 years, on the cuff - 1 year. The warranty does not apply to the handbag and packaging of the device. Extract from the Government Decree of 19.01.1988 No. 55. Approved "List of non-food goods of appropriate quality, not subject to return or exchange for a similar product of other size, shape, size, style, color or equipment." This list includes: "1 ... instruments, instruments and medical equipment ...". Information on maintenance, both under this warranty and paid, can be obtained from an authorized service-consulting point or by phone of a free hotline for Russia 8-800-200-33-22, or on the website of the company "Alfa-Medica" www.alpha-medica.ru. Service life of devices B.Well - not less than 10 years. The policy of the company "B.Well" provides for continuous improvement of products. In this regard, the company reserves the right to make full or partial changes to the products without prior notice and in accordance with production requirements.

12. DISCRIMINATION OF SYMBOLS USED FOR PACKAGING

SIGN OF COMPLIANCE

SIGN OF APPROVAL OF TYPE MEANS OF MEASUREMENTS

REFER TO THE INSTRUCTIONS FOR USE

IT IS POSSIBLE TO DISPOSE OF TOGETHER WITH HOUSEHOLD WASTE

MANUFACTURER

CE CERTIFICATION


MODEL

SERIAL NUMBER

APPLICATION: temperature limitation

STORAGE: storage temperature limit

Manufacturer: B. Well Swiss AG, Bahnhofstrasse 24, 9443 Widnau, Switzerland

Place of Manufacturer:

Wenzhou Bokang Instruments Co., Ltd., No. 1500 Haining Road, Haibin, Longwan, 325024, Wenzhou, China

Wuxi Exanovo Medical Instrument Co., Ltd. No. 42 Xixin Road, Zhangjing, Xibei Town, Wuxi City, 214194 Jiangsu, China.

The official importer, the organization authorized by the manufacturer, accepting the claims: Alfa-Medica JSC, 125493, Moscow, Avangardnaya str. 3, unit I, 4th floor, office 2402, phone: +7 (495) 645-86-99, e-mail: alphamed@online.ru

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
- | | |
|---|---------------------------------|
| 2. has been signed by
подписан | Thomas Gschwend |
| 3. acting in the capacity of
в качестве чего | legalizing officer |
| 4. Bears the seal/stamp of
Со штампом | Öffentliche Beglaubigungsperson |

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