Baxter Healthcare SA 8010 Zürich, Switzerland


Tel/Fax: +41 44 878 62 47 / +41 44 878 63 50

e-mail: steven_martin@baxter.com

Name: Steven Martin

Position: VP Quality EMEA

Signature:

Date:

Baxter Healthcare SA
Postfach
8010 Zurich

Instructions for use

OptiCap Disconnect Cap with Povidone-Iodine Solution

Description

This device contains two plastic disconnect caps for Automated Peritoneal Dialysis (APD), each containing povidone-iodine. The larger (FlexiCap) disconnect cap is in-tended to protect the Luer patient line connector of Baxter APD disposable sets. The smaller (MiniCap) disconnect cap is in-tended to protect the female Luer connector of the Baxter transfer set. Caps are intended for use only during dwell phase of therapy.

Contraindications

Do not use this product if there is a known history of allergic reaction to iodine. Contact your clinician for further instructions.

Warnings

Use aseptic technique. Contamination of any portion of the fluid path may result in peritonitis.

Do not use if pouch has been previously opened, damaged or stained on the outside of the pouch, or sponge is dry. Discard if sponge is touched.

Store at room temperature. Avoid excessive heat.

This product is intended for single use only. Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.

Keep out of the reach of children. In case of accidental ingestion, contact a physician immediately.

It is recommended that thyroid function be monitored in patients with small peritoneal dialysate fill volumes, typically infants and children. In order to minimize iodine expo-sure, drain the contents of the peritoneal cavity to the drain option prior to initiation of the next fill cycle whenever clinically possible.

Many dialysis products available from other manufacturers are used with equipment or disposables from Baxter Healthcare Corporation. Baxter has no control over variability, tolerances, mechanical strength or changes in these products which may be made from time to time. Therefore, Baxter cannot ensure that the dialysis products of other manufacturers, when connected with its products, will function in a satisfactory manner.

Directions for Use

Disconnection during Dwell

1. Close the transfer set.
2. Close the clamp on the patient line of the disposable set.

Caution: Use aseptic technique.


3. Follow the training provided by your clinician. Wear a face mask. Wash and dry your hands or use a disinfectant. Dry hands completely.
4. Open the **OptiCap** disconnect cap pouch.
5. Disconnect the transfer set from the patient line of the disposable set.
6. Place the patient line in the organizer.
7. Hold the transfer set in a downward position and immediately attach the **MiniCap** disconnect cap.
8. Tighten the **MiniCap** disconnect cap until fully secured.
Note: Avoid over tightening the **MiniCap** disconnect cap.
9. Attach the **FlexiCap** disconnect cap to the patient line connector in the organizer.
10. Tighten the **FlexiCap** disconnect cap until fully secured.

Reconnection and Continuation of Therapy

Caution: Use aseptic technique.

1. Follow the training provided by your clinician. Wear a face mask. Wash and dry your hands or use a disinfectant. Dry hands completely.
2. Remove the patient line of the disposable set from the organizer.
3. Remove the **FlexiCap** disconnect cap from the patient line connector of the disposable set.
Note: If povidone-iodine is not visible on the patient line connector, discard disposable set.
4. Remove the **MiniCap** disconnect cap from the transfer set connector.
5. Immediately connect the patient line connector of the disposable set to the patient transfer set.
6. Open the patient line clamp.
7. Open the transfer set.
8. Discard both disconnect caps in accordance with local waste disposal guidelines.

Design of the medical product

Application area

A disposable cap product used during disconnection situations of peritoneal dialysis therapy.

Toxicological

Material contained in key functional arcsos	
Component	Material
FlexiCap	Low density polyethylene (LDPE), white
MiniCap	Low density polyethylene (LDPE), white
Spacer	Polypropylene, Light blue
Sponge	Polyurethane foam, white
Pouch	Laminated Aluminum Foil Pouch
Mailer	Kraft mailer carton

Symbols

Symbol	Transcript	Symbol	Transcript
	The Manufacturer: Baxter Healthcare SA 8010 Zurich, Switzerland		Do not re-use
	Use-by date		Batch number
	European Conformity		Contains Povidone-Iodine
	Consult instructions for use		Keep away from heat
 RT	Store at room temperature		Do not use if package is damaged

Manufacturer's warranty

Manufacturer guarantees the conformity of medical device compliance with shipping, storage and use transportation to medical device

Date of Expiry

Subject to the rules of storage and transportation of the expected shelf life is one year (12 months).

Address of the manufacturer & addresses of production sites:

Legal manufacturer:

Baxter Healthcare SA
8010 Zurich
Switzerland

Manufacturing site:

Baxter Healthcare S.A.
Foxford Road, Swinford, Co.Mayo,
Ireland

Address of the representative in the territory of the Russian Federation:

125171, Moscow,
Leningradskoe shosse, 16A, building 1.

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