I, Michelle Metzger, hereby swear (or affirm) that the attached **Appendix to the Instruction for use** is a true copy of the original document.

Michelle Metzger

3.26.15

Michelle Metzger (ROS Name)
3850 Brickway Blvd. BWB2
Santa Rosa, CA 95403

Date

Country: Russia

Re: Piton GC Carotid Guide Catheter

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California

ss.

County of Sonoma

Subscribed and sworn to (or affirmed) before me on this March 26th, 2015, by Michelle Metzger proved to me on the basis of satisfactory evidence to be the person who appeared before me.

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Information about sticker

The final packaging for the Russian market is accompanied with additional sticker. Devices supplied in packaging marked with the following possible information:

- name of medical device in Russian language;
- name of the manufacturer;
- manufacturing sitecode UPN (unique product number);
- code CFN (customer facing number);
- date and number of Registration Certificate;
- sign of conformity in Russia ;
- warning and informing symbols and inscriptions.

Sample of the sticker

	Product name	
	Name of the manufacturer	Code UPN
	manufacturing site	Code CFN
	See IFU, Date and Number of RC	

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STERILIZATION

The catheter is sterilized in Ethylene Oxide and will remain sterile provided the packaging is not opened or damaged.
 For single use only. Do not re-use.

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DISPOSAL

After use, this products may be a potential

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
I, Beibei Sun, hereby swear (or affirm) that the attached Instructions for Use is a true copy of the original document.

Beibei Sun (RAS name)
3850 Brickway Blvd. BWB2
Santa Rosa, CA 95403

Date Aug 12, 2014


Russia / Piton GC


State of California

ss.

County of Sonoma

Subscribed and sworn to (or affirmed) before me on this August 12th, 2014, by Beibei Sun proved to me on the basis of satisfactory evidence to be the person who appeared before me.

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





DEVICE NAME

Invatec Piton™ GC is a steerable carotid guide catheter manufactured by Invatec S.p.A.

DESCRIPTION

Invatec Piton™ GC Carotid Guide Catheter is a dedicated guide catheter for supraaortic application. Proximally the shaft is reinforced by braiding. The pre-shaped distal part is reinforced by a flat-coil and ends in a radiopaque, soft tip. The short transition part between braided shaft and coiled tip incorporates the side hole for insertion of a second, stiff guidewire. Please refer to the single product label for precise product specifications.

INDICATIONS

The Piton™ GC Carotid Guide Catheter is intended to facilitate the introduction and placement of interventional devices (e.g. guidewires, stent delivery systems, dilation balloons, angiographic- or micro-catheters, etc.) into the human vasculature to treat vascular obstructive disease, including but not limited to the supraaortic vessels.

CONTRAINDICATIONS

There are no known contraindications for this device.

WARNINGS

- This device is designed and intended for single use only. DO NOT RESTERILIZE AND/OR REUSE. Reuse or resterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including (but not limited to) the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Reuse or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness and death. INVATEC will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- Inspect the device, prior to procedure, to verify functionality and lack of damaged parts. Do not use the device if the outer or the inner package is damaged or opened.
- Do not expose the device to organic solvents, e.g. alcohol.
- When the guide catheter is in the body, it should be manipulated while under sufficient and/or high quality fluoroscopy. Do not attempt to advance the catheter without paying careful visual attention to the tip response.
- Only physicians thoroughly trained and educated in the performance of transluminal angioplasties should use this device. Physicians should keep themselves informed and updated on recent publications about transluminal angioplasty techniques.
- Do not use with Lipiodol contrast media, or other such contrast media, which incorporate the components of these agents.
- If resistance is felt during insertion and/or manipulation of the device, determine the cause of resistance under fluoroscopy and take immediate remedial action.
- If the guide catheter is moved excessively, it may be damaged or may dislodge debris or emboli from the vessel wall, which can result in embolization of distal vessels.
- When torquing the guide catheter inside the vessel, do not exceed two rotations (720°) in the same direction without movement of the tip. This may damage the guide catheter or separate the tip, resulting in vessel injury or fragments might embolize distal vessels.
- Use the device prior to the "Use By" date specified on the package.

PRECAUTIONS

- Use prior to "Use By" date.
- Appropriate drug therapy (anticoagulant, vasodilator, etc.) should be administered to the patient according to standard protocols for percutaneous transluminal angioplasty.
- Exercise care during handling in order to avoid possible damage to the guide catheter. Avoid acute bending or kinking of the guide catheter itself. Do not use a guide catheter that has been damaged. Do not attempt to straighten a guide catheter that has been kinked or bent. Do not advance any device through a kinked guide catheter to reduce the potential of entrapment or device damage.
- Keep at controlled room temperature, in a dry place. Keep away from sunlight.
- Should any resistance be felt at any time during manipulation, insertion or withdrawal of the guide catheter, do not force or continue: stop the procedure immediately and determine the cause of the resistance before proceeding. Excessive force against resistance may result in damage of the device or the vessels.
- Various devices and catheters are designed to be compatible with different size guide catheters. Carefully check and match the device compatibility to the proper inner diameter of the guide catheter prior to use.
- Only interventionalists who have sufficient experience should carry out transluminal angioplasty. A thorough understanding of the technical principles, clinical applications and risks associated with transluminal angioplasty is necessary before using this device.
- Allergic reactions to contrast medium should be identified before treatment.
- Check that the catheter size is compatible with the interventional devices (e.g. guidewires, stent delivery systems, dilation balloons, angiographic- or micro-catheters, etc.) for the procedure which are intended to be used in.

POTENTIAL COMPLICATIONS / ADVERSE EFFECTS


Complications associated with the use of a Piton™ GC Carotid Guide Catheter are similar to the ones associated with standard PTA procedures. Possible complications may include, but are not limited to:

- Death.
- Stroke.
- Vascular spasm.
- Vessel dissections, perforation, rupture.
- Hypo- / hypertension.
- Sepsis / Infection.
- Arrhythmia, including ventricular fibrillation.
- Hemorrhage or hematoma.
- Arteriovenous fistula.
- Pseudo-aneurysm.
- Thrombosis and / or embolization.

DIRECTIONS FOR USE

Carefully remove the guide catheter from its packaging.

1. Prior to procedure, carefully examine all equipment to be used, including the guide catheter, to verify proper function. Check that the catheter and sterile packaging have not been damaged during

- shipment and that the catheter size is compatible for the specific procedure for which it is intended.
2. Flush the guide catheter with heparinized saline.

INSERTION TECHNIQUE

1. Prepare the patient to angioplasty using standard technique. Appropriate anticoagulation and vasodilation therapy should be used.
2. Under angiographic control, advance a stiff 0.035" exchange guidewire (300 cm long, called pivotal wire) through an at least 8F introducer sheath into the aortic arch.
3. Insert the proximal end of the pivotal wire into the side-hole of the Piton™ GC Carotid Guide Catheter (GC).
4. Before introducing the Piton™ GC through the introducer sheath, straighten the tip of the Piton™ GC by inserting a 0.035" guidewire of choice for Carotid access (e.g. floppy, hydrophilic guidewire, called standard wire) from the proximal hub of the GC.
5. As soon as the tip of the Piton™ GC is straightened, insert the standard wire into the introducer sheath as well and follow with the Piton™ GC together with the standard wire in place.
6. Under angiographic control advance the entire system over the pivotal wire into the aortic arch until the side hole of the Piton™ GC is in front of the ostium (e.g. Common Carotid Artery) of the vessel to be treated. The standard wire has to exceed the Piton™ GC tip as well as needed in order to always keep it straightened.
NOTE: the position of the side hole is located between the coiled and the braided part of the Piton™ GC tip, which can be angiographically identified as a radiopaque gap between the two structures.
7. Pull back the standard wire partially, so that the tip curve of the Piton™ GC will start to reshape.
8. When the tip curve has been reshaped, take the Piton™ GC at the hub and rotate it carefully in order to engage the tip of the Piton™ GC into the ostium of the vessel.
9. When the tip has been engaged into the ostium, carefully re-advance the standard wire to an appropriate level distal of the ostium.
10. Pull back the pivotal wire just proximal of the side hole, in order to leave the external position of the Piton™ GC. Re-advance the pivotal wire through the inner lumen of the Piton™ GC up to the same level of the standard wire.
11. Under angiographic control advance the Piton™ GC over the two guidewires (pivotal and standard wire) into the vessel, until it reaches a stable position.
12. Retract the standard wire and start the planned intervention.

HOW SUPPLIED

The Piton™ GC Carotid Guide Catheter is supplied sterile and intended for single use only. It is sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. Use the product latest by the printed "Use by" date.

CAUTION: Do not use if the inner package is open or damaged.

STORAGE

Store at controlled room temperature in dry place. Keep away from sunlight. Do not expose to organic solvents (e.g. alcohol), ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used prior to the expiration date on package label.

WARRANTY / LIABILITY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in Invatec's instructions for use are expressly to be considered as an integral part of this provision. Invatec guarantees the product until the expiry date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. Invatec disclaims any warranty of merchantability or fitness for a particular purpose of the product. Invatec is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on Invatec's part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of Invatec, however originating, from the product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Invatec's legitimate interest in limiting its liability or warranty. No person has any authority to bind Invatec to any warranty or liability regarding the product.

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