«Approved»

MEDIVATORS
A Cantel Medical Company

MEDIVATORS Inc.
14605 28th Avenue North
Minneapolis, MN 55447 U.S.A.

Call: (800) 327-2340
Fax: (800) 660-8498

Megan Dickey

Signature \ Seal

INSTRUCTION FOR USERS

«Endoscopy CO2 insufflator Endo Stratus CO2 Insufflator»

1. Product's description

Endoscopy CO2 insufflator Endo Stratus CO2 Insufflator consist of:

- Endoscopy CO2 insufflator Endo Stratus CO2 Insufflator EGA-501E
- Power cord
- IFU

2. Indication for use

The ENDO STRATUS CO2 insufflator is designed to use CO2 as a distention media in the gastrointestinal tract when used in conjunction with a gastrointestinal endoscope during endoscopy diagnostic and treatment and for maintaining the nominal pressure.

The ENDO STRATUS™ CO2 insufflator contains a CO2 warming feature that allows the clinician to modulate and control delivery of CO2 to the patient that has been warmed to a target temperature 37°C.

3. Definitions:

- Throughout this document the ENDO STRATUS™ CO2 insufflator may be referred to as “unit” or “device”
- L/min – Liters per minute (flow)
- kPa – kilo Pascal (pressure)
- VAC – Volts Alternating Current (electrical potential)
- W – Watts
- Hz – Hertz (frequency)
- Mm – Millimeters
- Kg - Kilogram
- A/m – Ampere per meter

Signal words:

WARNING:

Indicates a potentially hazardous situation which if not avoided, could result in death or serious injury.

CAUTION:

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.

4. Set:

Endoscopy CO2 insufflator Endo Stratus CO2 Insufflator consist of:

- Endoscopy CO2 insufflator Endo Stratus CO2 Insufflator EGA-501E, 1pc
- Power cord, 1pc
- IFU, 1pc

5. Description of the product

Theory of operation

The ENDO STRATUS CO2 insufflator works by regulating CO2 from a tank source to a preset pressure and then controlling the output of CO2 to an GI endoscope system for ultimate delivery as a distention media in the gastrointestinal tract. The clinician will use the air/water valve on the GI endoscope and visual feedback of the GI endoscopic system to manually distend the gastrointestinal tract with CO2.

The ENDO STRATUS CO2 insufflator contains a CO2 warming feature that allows the clinician to modulate and control delivery of CO2 to the patient that has been warmed to a target temperature of 37°C with an upper tolerance of +3°C. Warming is achieved as CO2 flows across a heating element which can be switched on or off by the user,

depending on if warming of the CO₂ is desired. The heating of the CO₂ is accomplished through an internal 20W heating element with dual redundant temperature sensors.

The flow rate can be adjusted to low, medium or high settings depending on the desired rate of insufflation, with flow rates of 1.4 L/min, 2.4 L/min and 3.5 L/min maximum, respectively.

When the unit is turned on, the CO₂ flow and heating function are off by default. To initiate gas flow the user must turn the flow setting knob from the “off” position. The unit will detect gas flow and will not allow the heater to operate unless gas flow is detected. In addition, there is a low gas pressure sensor and the unit will indicate when the input CO₂ pressure reaches 172,4 kPa by illuminating a yellow light on the front panel.

The unit also has an output to power the water bottle warmer. The water bottle heating element maintains the temperature of the water in a sterile water bottle at 37±3°C. It is controlled by dual temperature sensors for redundancy and safety.

The heater of water bottle used for support constant temperature at 37±3°C at water bottle. To connect water bottle through consumables to endoscope and irrigation pump, through “water” button on an endoscopy, water from water bottle goes to endoscopy, to make an irrigation of a necessary area. Heater function is an additional function and physician could use any water bottles that are approved to use on medical centers in Russian Federation. Endoscopy CO₂ insufflator Endo Stratus CO₂ Insufflator doesn't consist of a water bottle. Water bottle holder holds a water bottle.

6. Application area

Insufflator is used for endoscopic treatment and biopsy in endoscopic department

7. Contradictions

The ENDO STRATUS™ CO₂ insufflator should only be used by or under the direct guidance of a physician experienced in the standard practices of gastrointestinal endoscopy procedures. It should only be used for procedures where insufflation of the gastrointestinal tract is necessary to support navigation of the GI endoscope and performance of any necessary evaluation procedures. The device should not be used in any other method or for any other treatments or procedures.

The device is contraindicated for laparoscopic or hysteroscopic insufflation. It must not be used for intrauterine distension.

The device is contraindicated for CT colonography. Also, the device is not intended for use with or exposure to magnetic resonance imaging (MRI) systems. Do not use the device in an environment where it could be exposed to MRI systems.

8. Possible side effects

Possible side effects of using insufflator associated with possible side effects of endoscopic procedures.


9. Service regulations

Insufflator should be used at endoscopic department of the hospitals.


Initial setup

1. Place the ENDO STRATUS CO₂ insufflator on a flat surface such as an accessory cart or other suitable working surface.
2. Before beginning, inspect the ENDO STRATUS CO₂ insufflator for any signs of damage.
3. Connect the power cord to the back of the device.
4. Before connecting the power cord to the hospital grade wall receptacle, make sure the power is off, and no accessories are connected. Then plug the power cord into the wall receptacle.

5. Install the pump head onto the device by aligning the flat shaft on the front of the device with the slotted hole in the pump head and rotate the pump head until it locks in place.
6. Assemble the water bottle holder and warmer assembly in the following manner:
 - A. The following supplies are necessary and can be found in a bag attached to the water bottle holder:
 - a) Short Screws (2)
 - b) Long Screws (2)
 - c) Large Hex Key
 - d) Small Hex Key

- b. Attach the water bottle holder to the base tray using the (2) small screws and small hex key.

- c. Place the ENDO STRATUS™ CO2 insufflator on the base tray and install the (2) long screws using the large hex key.

Adjust the water bottle hold-down bracket to fit the water bottles used by your facility. This is done by loosening the nut on the back and sliding the bracket up or down so it contacts the water bottle and then re-tightening the nut.

Plug the water warmer cable into the receptacle marked "water warmer" on the back of the unit.

Air hose attachment

Tools required: (not included)

- 9/16" open-end wrench
- 3/4" open-end wrench

Accessories required:

- "D" or "E" size CO2 tank filled with medical grade CO2 (not included)
- High pressure hose

1. Assemble the high pressure hose and the yoke by tightening the hose with the 3/4" wrench while holding the yoke with the 9/16" open-end wrench.
2. Unscrew the T-handle on the yoke to allow it to slip over the top post of the CO2 tank. The correct orientation is so that the two pins on the yoke line up with the two alignment holes on the tank. (CO2 tank is not included with system, purchase separately).
3. After aligning the two pins, tighten the T-handle firmly.
4. Tighten the other end of the high pressure hose and yoke assembly to the back of the ENDO STRATUS™ CO2 insufflator. Hold the brass fitting with wrench and use a second wrench to tighten the hose nut. **DO NOT USE TEFLON TAPE OR THREAD SEALANT ON THE FITTINGS.**

Pre-procedure setup

1. Open the valve on the CO2 tank approximately 1 turn.
2. Turn the power switch to "On".
3. Test the setup by turning the flow control knob position clockwise. Verify that there is CO2 flowing out the front output of the unit. If there is no CO2 observed out the front of the unit verify that the CO2 tank is full and check the connection outlined in air hose attachment above.
4. Turn off the flow control knob.
5. Connect the tubing to the front output of the device. Use only ENDO SMARTCAP® tubing or ENDOGATOR® hybrid tubing manufactured by MEDIVATORS.
6. When using ENDO SMARTCAP® CO2 source tubing, attach ENDO SMARTCAP® source tubing to ENDO SMARTCAP® tubing or ENDOGATOR® hybrid tube set. Attach ENDO SMARTCAP® tubing or ENDOGATOR® hybrid to the GI endoscope according to the manufacturer's instructions.

Water Pre-warming

1. The ENDO STRATUS CO2 insufflator can be used to control the water bottle warmer system. The warmer system can maintain the temperature of bottles of sterile water at 37°C (±3°C). However, the system is not intended to raise the temperature of the water; therefore, pre-warming of the water bottle is required.

CAUTION:

The system is not indicated for warming of a water, pre-warming is required!

2. The water bottle must be pre-warmed with a method such as a fluid warmer if warm water is desired. Water bottle pre-warming should be performed according to the following instructions:

Place sterile water bottle in a fluid warming oven set to 37°±3°C. Do not exceed this temperature. If warm water is not desired, turn the water warmer on the pump off and place a bottle of room temperature sterile water on the water bottle holder.

NOTE:

If the water warmer has been enabled and has exceeded the upper temperature set point limit, the water warmer's switch indicator light will turn blue and the warmer will automatically shut off. If the switch indicator light is flashing blue, ensure the water warmer cable is plugged into the back of the unit. If it is connected and the light is flashing, this may indicate a problem with the unit. Contact customer service at 800-444-4729.

Read and become familiar with all manufacturer's instructions on warming ovens and water bottles regarding maximum fluid temperatures. Never exceed 40°C during water pre-warming. Never use a microwave oven to warm a water bottle, as this could heat the water to a dangerously high temperature or could result in uneven warming.

CAUTION:

The surface of the water bottle warmer may be hot to the touch.

Operation

1. After the pre-procedure setup has been completed as outlined in section VIII, flow can be initiated.
2. Turn the flow control knob to the first position. "L" will illuminate for "LOW". If a higher flow rate is desired, turn the knob to the next position for "MEDIUM" and one more position for "HIGH". "M" will illuminate for "MEDIUM" and "H" will illuminate for "HIGH". The knob will turn freely in each direction, and will increase or decrease the flow rate accordingly. See section V for flow rates.
3. If warming of the CO₂ is desired, press the button for the CO₂ heater. The button will illuminate green.
4. If the device is attached to a water bottle warmer, and water warming is desired, press the button for water bottle warmer. The button will illuminate green.
5. Use the appropriate air/water valve on the GI endoscope and visual feedback of the GI endoscopic system to manually distend the gastrointestinal tract with CO₂.

Shutdown

1. When the procedure is over, turn the knob counter clockwise until all lights are off, turn off the water bottle warmer (if attached) and turn off the CO₂ heater.
2. Turn off the main power switch and disconnect the tubing from the front of the unit.
3. Discard the tubing according to the manufacturer's instructions for proper disposal.

Front panel controls/connectors**NOTE:**

Note: All status indicator lights will be green during normal operation.

1. **On/Off push button**
Turns main power on or off to the unit
2. **Water warmer On/Off**
Push to turn water warmer on or off
3. **CO₂ Flow control**

Controls flow rate of CO₂ output. The three levels are represented by three green indicators of "L", "M" and "H" which are illuminated to indicate the flow rate setting.

The flow rates are as follows*:

Low – 1.4L/min

Medium – 2.4L/min

High – 3.5L/min


**flow rates shown are maximum free flow measurements at output of unit*

Low gas pressure is indicated by the words "LOW PRESSURE" flashing to notify the user that the gas input pressure has fallen below 72,4 kPa.

4. **CO₂ gas output connector**
Connect to tubing and GI endoscope.
Standard male luer lock connector.

5. CO2 gas heater On/Off

Push to toggle CO2 gas heater On or Off

(the appearance of your device may differ slightly from picture shown above)

Rear panel connections

1. Water bottle warmer connection

Power connection for water bottle warmer

Refer to section "service regulations" for warmer explanation

2. CO2 input connection

For connecting a "D" or "E" sized CO2 tank. Use only pre-filtered medical grade CO2 gas.

3. Equipotentiality (Earth Ground)

4. Fuse receptacle

WARNING:

Replace only with type and rating marked M10AL250V

After fuse receptacles were used, you should change it to fuse receptacles with the same characteristics (M10AL250V). You should push the holder of fuse receptacles, see figure 1.

Figure 1. Holder Fuse receptacle (open)

After that you should change fuse receptacles (2 pieces), you should push on it and change to a new one and close the holder (see figure 2)

Figure 2. Holder Fuse receptacle (close)

5. AC power connection

AC power input 100-240 VAC, 50-60HZ


WARNING:

Use only the hospital grade power line cord supplied with this unit. Connect only to a power receptacle marked as hospital grade.

6. Water bottle holder

The water bottle holder contains the warming system for the water bottle and connects to the main ENDO STRATUS CO2 insufflator by means of an electrical cable which plugs into the receptacle marked "water warmer".

(the appearance of your device may differ slightly from picture shown above)

10. Classification of the medical device

Endo Stratus CO2 Insufflator is a non-invasive medical device used for insufflations of CO2 in in the gastrointestinal tract. The insufflator is a non-sterile, reusable device.

11. Safety measures

WARNINGS

- A. To reduce risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel.
- B. To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth..

- C. The ENDO STRATUS CO2 insufflator is not suitable for use in the presence of a flammable anesthetic mixture with oxygen..
- D. The ENDO STRATUS CO2 insufflator shall be sold only by prescription for use by physicians/clinicians who are trained and qualified.
- E. Only qualified medical personnel in an acceptable medical facility should operate the ENDO STRATUS CO2 insufflator.
- F. The ENDO STRATUS CO2 insufflator should be connected to a properly grounded receptacle otherwise grounding reliability cannot be achieved.
- G. Extreme precaution must be taken when handling liquids around electrical equipment. DO NOT operate the ENDO STRATUS™ CO2 insufflator if liquid has been spilled on the unit.
- H. Never place or stack the ENDO STRATUS™ CO2 insufflator on other electrical equipment other than the ENDO STRATUS™ irrigation pump. Electromagnetic or other interference may occur between the ENDO STRATUS™ CO2 insufflator and other electronic devices. The equipment or system should be observed to verify normal operation in the configuration in which it will be used.
- I. The ENDO STRATUS™ CO2 insufflator should only be used in conjunction with other equipment whose safety against leakage currents has been established.
- J. The instructions contained in the operating manuals of any equipment to be used in conjunction with the ENDO STRATUS™ CO2 insufflator must be followed to avoid any possible hazard from incompatibility.
- K. The instructions for use described in this manual MUST be followed. Otherwise, compromised safety, malfunction, injury to the operator and/or patient, or costly damage to the unit and other equipment may occur.
- L. The ENDO STRATUS™ CO2 insufflator must be connected to an appropriate power source when loss of power source would result in an unacceptable risk.

CAUTION

- A. If emergency or abnormal function occurs, immediately turn off power to the unit.
- B. Read and understand all warnings that come with your commercially available “D” or “E” sized CO2 tanks.
- C. Remove power from the device before initiating any field servicing of the replacement parts.
- D. There are no user serviceable parts inside this unit. Repairs to the ENDO STRATUS™ CO2 insufflator should only be performed by qualified service personnel.
- E. When the low pressure warning indicator illuminates, exchange the CO2 tank as soon as possible to avoid loss of function.
- F. Always keep a spare full tank of CO2 nearby for quick access.
- G. Always keep the CO2 tank in an upright position to avoid fluid entering the unit.
- H. Do not use the device if the enclosure is damaged or enclosure integrity has been compromised.
- I. Do not attempt to operate the device before reading and understanding all sections of this manual.
- J. Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in section “safety measures”.
- K. Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment. Do not expose the device to sources of electromagnetic interference such as CT equipment, diathermy equipment, cellular phones, RFID tags and metal detectors.

12. Technical specifications

Electrical specifications

Input Voltage:: 100-240 VAC

Input Frequency:: 50-60 Hz

Power Consumption:: 82 VA

Fuse Rating:: M10AL250V

Medium acting, 10 amp, low breaking capacity, 250 volt

Replace fuses only with those of same type and rating

Certifications: IEC-60601-1, IEC-60601-1-2, IEC-60601-2-18

Classification: Class 1 Type BF

IP rating (Ingress Protection): IP24

WARNING:

Grounding reliability can only be achieved when connected to a receptacle marked "Hospital Grade"

Mechanical characteristics

Physical Dimensions::

Height:	121 mm
Width:	197 mm
Depth:	349 mm
Weight:	4,8 kg

CO2 Specifications

Free Flow Rates (maximum unrestricted output to GI endoscope system):

Front panel setting

Low:	1,4 L/min
Medium:	2,4 L/min
High:	3,5 L/min

GI Endoscopic Flow Rates (approximate flow rate at GI endoscope output):

Front panel setting

Low:	0,6-0,9 L/min*
Medium:	1,2-1,4 L/min*
High:	1,5-1,8 L/min*

The difference between free flow rate and flow rate that was set (gas consumption):

Front panel setting

Low:	0,5-0,8 L/min
Medium:	1,0-1,2 L/min
High:	1,7-2,0 L/min

*These values are approximate and are based on testing with respective models of GI endoscopes. User results may vary depending on GI endoscope used, channel diameter, and channel length.

Output Connection:: Male luer lock

Input Connection:: 1/4" male flare fitting.

CO2 Heater:: 20W, internal, flow through, with dual redundant temperature sensors.

CO2 Output temperature:: 37⁰C (+3⁰C)

Operating Pressure::

Input Pressure: 13 100kPa± 10kPa (max)
172,4kPa± 0,5kPa (min)

**The low pressure warning light will illuminate at 25 PSI, although the unit may still flow. It is recommended to replenish the CO2 supply when the light illuminates.*

Pressure regulator set point: 55,16kPa (min)

Internal safety relief valve: ≤ 82,74kPa

Output pressure ratio: 75 kPa to 106 kPa

Power of water heater: not more than 20 W

Environmental requirements:

Operating Temperature::	16 ⁰ to 24 ⁰ C (or +61 ⁰ F to 75 ⁰ F)
Operating Relative Humidity::	30% to 70% non-condensing
Operating Pressure::	75 kPa to 106 kPa
Storage Temperature::	-50 ⁰ to 50 ⁰ C (or -58 ⁰ F to 122 ⁰ F)
Storage Relative Humidity::	10% to 98% non-condensing
Storage Pressure::	70 kPa to 106 kPa

Products features:

- Internal CO2 warmer eliminates cold CO2.
- Includes control for optional water bottle warmer to maintain water at close to body temperature.
- Stackable with the ENDO STRATUS irrigation pump
- Standard male luer lock output
- Uses standard "D" or "E" size CO2 cylinders
- Accepts input from wall CO2 with standard DISS type connection
- Usable at CO2 input pressures lower than most comparable units, making efficient use of tank supply
- Works with disposable ENDO SMARTCAP tubing and ENDOGATOR hybrid tubing for safety and compliance

Safety features:

- Dual inline pressure regulators ensure a constant pressure
- Additional mechanical pressure relief at $\leq 82,75$ kPa prevents overpressure
- Low CO2 input pressure is indicated by a yellow light on the front panel
- On initial power up, gas flow is off until the flow knob is turned on.
- CO2 warmer will not operate unless flow is initiated.

Materials and marks used for manufacturing of housing of the device and power cord:

№	Part	Materials	Trademark	Manufacturer
1	Housing of the device	ABS plastic	ULTEM 1010-7101	«Polimerland Inc», USA
		Stainless steel	UNS S44100	«AK Steel», USA
2	Power cord	Polypropylene	PRONAT.03X24X47L/S	ePlastics Inc., USA
		Brass	223 Brass cleaned	USA Brass Company, Inc., USA

13. Service

The device is not user serviceable (excluding the replacement parts listed in section V). Please contact Customer Service and obtain a Return Material Authorization (RMA) Number. The RMA number must accompany all units and paperwork for proper service or warranty work.

14. Care and maintenance

Cleaning

- Before cleaning, ensure the power to the unit is turned off and the electrical cord unplugged.
- The outside surface of the unit may be cleaned with a damp cloth, 70% isopropyl alcohol solution, or a 10% bleach with water solution as often as deemed necessary.
- To disinfect the outside surface, use a mild disinfectant as needed according to the manufacturer's instructions..
- Do not use abrasive or harsh cleaners..
- Do not allow liquid to enter the unit.
- Do not sterilize the unit by any means.

Technical support

Materials and workmanship of this product, ENDO STRATUS CO2 insufflator will be warranted for a period of one (1) year from the original date of purchase. If this product should become inoperable due to a defect in materials or workmanship during this one (1) year warranty period, the product will be repaired or replaced.

Change of fuse receptacles

After fuse receptacles were used, you should change it to fuse receptacles with the same characteristics (M10AL250V). You should push the holder of fuse receptacles, see figure 1.

Figure 1. Holder Fuse receptacle (open)

After that you should change fuse receptacles (2 pieces), you should push on it and change to a new one and close the holder (see figure 2)

Figure 2. Holder Fuse receptacle (close)

15. Packaging

A single package consist of: cardboard box with cardboard cards and foam filler.

16. Marking

A single package marked by:

- The manufacturer,
- Address and phone number of the manufacturer,
- Description of the package
- LOT
- Storage temperature
- Reusable device
- Avoid a direct sun rays

-Distributor

-Date of manufacturing

- Warranty

- Fragile

An additional information could be added


- follow the IFU

- Relative humidity

- Pressure

Symbols and warnings, mentioned on the device and label:

Symbols and warnings	Description
 / model	Model
	Serial number
 Manufactured by	Manufacturer
Water Warmer 	Please, use an instruction (water warmer)
Danger. Equipment not suitable for use in the presence of flammable anesthetic mixture with oxygen (Danger: l'équipement n'est pas conçu pour être utilisé en présence d'anesthésiques inflammables enrichis en oxygène)	Danger. Equipment not suitable for use in the presence of flammable anesthetic mixture with oxygen
CO2 Input 	Please, use an instruction (CO2)
 1900 PSI/131 Bar Maximum. Use medical grade Co2 gas only (Maximum 1900 PSI/131 Bar Utiliser un gaz CO2 de qualité médicale seulement)	Warning. 1900 PSI/131 Bar Maximum. Use medical grade Co2 gas only
 Class 1 Type B IP24	Class 1 Type B defender class of the housing is 24
	Equipotentiality
 Warning: Grounding reliability can only be achieved when connected to a receptacle marked «Hospital Grade». (Avertissement: la mise à la terre	Warning: Grounding reliability can only be achieved when connected to a receptacle marked «Hospital Grade».

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17. Storage and transportation

Storage

Please, storage in a dry place and avoid high temperature, humidity, ventilation, sun rays, dust and air with an salt etc.
Please, check a possible defects of the package after a transportation.

Avoid places with chemicals and gases.

Avoid wet places.

The insufflator is acceptable for the temperature and humidity of the air during the transportation and storage according to group 5 (OG4) GOST 15150.

Transportation

The insufflator is placed to appropriate package with IFU. Transportation by all types of covered transport is allowed according to temperature mode of group 5 (OG4) GOST 15150.

18. Recycling and environmental requirements

All medical devices, packaging materials and accessories should be utilized according to state standards and rules.

19. Operations conditions

The Insufflator should be used at endoscopic department of the hospitals.

20. List of international and national standards

Meets standards international and Russian requirements: GOST RF 2.601-2013, GOST R 50444-92, GOST R 50267.0.2-2005, MU 287-113-98, GOST 15150, GOST R 50267.0-95, RD 50-707-91, GOST R 50267.18-94 Directive 93/42/EEC, ISO 13485:2012.

Meets the standard Directive 93/42/EEC dated 14.11.2013, appendix 2 to medical devices class 2a

21. List of international certificates and declarations of conformity

- Certificate ISO 13485:2003 dated 15.04.2014r.
- Directive EC 93/42/EEC dated 14.11.2013 r.
- Declaration of conformity dated 19.06.2014 r.
- Certificate to foreign governments dated 07.05.2014

22. Warranties

The warranty period is 1 (one) year. If this product should become inoperable due to a defect in materials or workmanship during this one (1) year warranty period, the product will be repaired or replaced.

The compliance with the requirements of technical documentation are guaranteed by the manufacturer if the conditions of using, transportation and storage were met.

For all questions, please, contact an official distributor in Russia – Endo Stars LLC

Legal address: 197101, Saint-Petersburg, Bolshaya Monetnaya str., 27, lit A, 12H

Tel/fax: 8-800-555-56-57

23. Table of EMC compatibility

Table of EMC compatibility

Table 1 – Guidance And MANUFACTURER'S Declaration ELECTROMAGNETIC EMISSIONS For All ME EQUIPMENT And ME SYSTEMS

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF Emissions CISPR 11	Group 1	Use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The insufflator are suitable for use in all establishments, other than domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes..
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	

24. Guidance And MANUFACTURER'S Declaration Electromagnetic IMMUNITY For All ME EQUIPMENT And ME SYSTEMS

Table 2 – Guidance And MANUFACTURER'S Declaration Electromagnetic IMMUNITY For All ME EQUIPMENT And ME SYSTEMS

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
Immunity Test	EN/IEC 60601 Test Level	Compliance level	Electromagnetic Environment – Guidance
ESD EN/IEC 61000-4-2	±6kV Contact ±8kV Air	±6kV Contact ±8kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%.
EFT EN/IEC 61000-4-4	±2kV Mains ±1kV I/Os	±2kV Mains ±1kV I/Os	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/ Dropout EN/IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	100% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles Note 1	Mains power quality should be that of a typical commercial or hospital environment. If the user of the insufflator requires continued operation during power mains interruptions, it is recommended that the insufflator be powered from an uninterruptible power supply or battery.
Power Frequency 50/60Hz Magnetic Field EN/IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Note 1 The EUT shuts off and must be manually restarted by the operator during a 5 second loss of AC Mains power.

25. Guidance And MANUFACTURER'S Declaration Electromagnetic IMMUNITY For ME EQUIPMENT And ME SYSTEMS That Are Not LIFE-SUPPORTING

Table 4 – Guidance And MANUFACTURER'S Declaration Electromagnetic IMMUNITY For ME EQUIPMENT And ME SYSTEMS That Are Not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity

Immunity Test	EN/IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	(3)Vrms	Portable and mobile communications equipment should be separated from the insufflator by no less than the distances calculated/listed below: $D=(3.5/V1)(\text{Sqrt } P)$ $D=(3.5/E1)(\text{Sqrt } P)$ 80 to 800 MHz $D=(7/E1)(\text{Sqrt } P)$ 800 MHz to 2.5 GHz Where P is the max power in watts and D is the recommended separation distance in meters.
Radiated RF EN/IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	(3)V/m	Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter..

26. Recommended Separation Distances Between Portable And Mobile RF Communications Equipment And The insufflators EQUIPMENT and SYSTEMS That Are Not LIFE-SUPPORTING

Table 6 – Recommended Separation Distances Between Portable And Mobile RF Communications Equipment And insufflators EQUIPMENT and SYSTEMS That Are Not LIFE-SUPPORTING

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $D=(3.5/V1)(\text{Sqrt } P)$	Separation (m) 80 to 800MHz $D=(3.5/E1)(\text{Sqrt } P)$	Separation (m) 800MHz to 2.5GHz $D=(7/E1)(\text{Sqrt } P)$
0,01	0,1166	0,1166	0,2333
0,1	0,3689	0,3689	0,7378
1	1,1666	1,1666	2,3333
10	3,6893	3,6893	7,3786
100	11,6666	11,6666	23,333

WARRANTY CERTIFICATE

Note of Sale

Product category _____

Manufacturer _____

Model _____

Serial/identification number _____

Date of sale _____

Seller organization name _____

Seller's signature _____

The product is tested, no defects found.

Operating rules were read and accepted by:

Buyer's signature _____

Warranty period – 1 years.

For warranty questions, please contact:

Official distributor in Russian Federation – Endo Stars LLC, legal address: 197101, Russia, Saint-Petersburg, Bolshaya Monetnay, 27, lit. A, 12H; tel./fax: +7-800-555-56-57, e-mail: endo@endo-stars.com.

ATTENTION! Make sure that the Warranty Certificate is filled in correctly. The serial/identification number must correspond exactly with the one in the Certificate. All fields in the Certificate (except the fields for the warranty center) must be filled in.

In case of inaccurate or incomplete filling of the Certificate, the warranty for the product is regarded as invalid.

The product should not bear the marks of mechanical or thermal action, as well as of various chemical agents. The product use is allowed only in accordance with the Operators Manual.

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