State of California
Secretary of State

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays / País:		United States of America	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Misti Brooks	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Deputy	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de		County of Los Angeles, State of California	
Certified Attesté / Certificado			
5. at à / en	Los Angeles, California	6. the le / el día	2nd day of April 2015
7. by par / por	Secretary of State, State of California		
8. N° sous n° bajo el número	49049		
9. Seal / stamp: Sceau / timbre: Sello / timbre:		10. Signature: Signature: Firma:	

This Apostille is the trilingual model Apostille Certificate as suggested by the Permanent Bureau and developed in response to the 2009 Special Commission on the practical operation of the Hague Apostille Convention.
This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.
This Apostille does not certify the content of the document for which it was issued.
This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.
To verify the issuance of this Apostille, see: www.sos.ca.gov/business/notary/apostille-search/.

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Los Angeles County Registrar-Recorder/County Clerk

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I, DEAN C. LOGAN, County Clerk of the County of Los Angeles, State of California, in and for said County DO HEREBY CERTIFY THAT **MIHO HIRAMATSU** at time of signing a duly commissioned, qualified and acting NOTARY PUBLIC, in the State of California, County of Los Angeles, empowered to act as such Notary in any part of this State and authorized to take the acknowledgement or proof of powers of attorney, mortgages, deeds, grant transfers, and other instruments of writing executed by any person, and to take depositions and affidavits and administer oaths and affirmations in all matters incident to the duties of the officer or to be used before any court, judge, officer, or board.

I FURTHER CERTIFY that the seal affixed or impressed on the attached document is the official seal of said Notary Public and it appears that the name subscribed thereon is the genuine signature of the person aforesaid, his (or her) signature being of record in this office.

(Valid only if the certification bears the embossed Seal of the Registrar-Recorder/County Clerk.)

IN WITNESS WHEREOF, I execute this certificate and have hereunto set my hand and affixed the seal of said County this

02nd day of April, 2015

DEAN C. LOGAN
Registrar-Recorder/County Clerk

By: Misti Brooks
MISTI BROOKS, Deputy County Clerk

BRF03-3/3/2009

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California

County of Los Angeles

On 12/1/2015 before me, Miho Hiramatsu
(Insert Name of Notary Public and Title)

personally appeared Takashi Suzuki

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature [Signature] (Seal)

SAKURA

VisionTek®

Digital Microscope Brightfield System

Operating
Manual

VisionTek®

**Operating
Manual**

Digital Microscope Brightfield System

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0006663-01 Rev. B

TABLE OF CONTENTS

Section	Page
1 INTRODUCTION	
Intended Use.....	1.1
General Description	1.1
Safety Precautions.....	1.1
Physical Characteristics.....	1.2
Cabinet Front (Figure 1-A).....	1.2
Delivery Contents	1.3
Cabinet Rear (Figure 1-C).....	1.4
PC Rear (Figure 1-D).....	1.5
Specifications.....	1.6
Models Covered By This Manual.....	1.6
Power Requirements	1.6
Power Ratings	1.6
Environmental Requirements	1.6
Dimensions.....	1.6
Instrument Life Expectancy	1.6
Safety Standards	1.6
2 INSTALLATION	
General Information	2.1
Environmental Factors.....	2.1
Unpacking.....	2.1
Unpacking and Installation.....	2.2
Delivery Contents	2.2
Work Surface.....	2.2
Warnings and Precautions.....	2.2
3 OPERATING INSTRUCTIONS	
Operating Instructions.....	3.1
System Start.....	3.1
System Operation	3.2
Software Operation for Single Live Slide Analysis.....	3.3
VisionTek Operation for Multi-slide Analysis.....	3.12
Digitized Image Files	3.15
Digital Recordings of Images.....	3.16

TABLE OF CONTENTS

	Viewing of Digital Slide	3.22
	Slide Tray – Slide and Data Management	3.24
	The Virtual Case Session	3.33
	The Digital Gallery	3.37
	Settings.....	3.43
	Import/Export Function.....	3.49
	Keyboard Shortcuts	3.50
4	RECOVERY PROCESS	
	Overview	4.1
	Actions for Recovering from an Error Condition.....	4.1
5	CARE OF THE INSTRUMENT	
	Maintenance As Needed.....	5.1
	Slide Tray cleaning	5.1
	Daily Maintenance	5.1
	Exterior Cleaning	5.1
	Annual Maintenance	5.1
6	TROUBLESHOOTING	
	General Troubleshooting Procedures	6.1
	Service Information	6.1
	When a Problem Occurs.....	6.1
	Where to Call for Service.....	6.1

INTRODUCTION

Intended Use

The VisionTek® Digital Microscope is used to visualize embedded or frozen human and animal tissue placed on a microscope slide. The image is viewed on a standard computer screen in real time with the option to store the field of view. The device is used in the same manner as light microscope when it is determined that a surgical specimen requires microscopic examination for pathology by a qualified pathologist. Like the light microscope, the VisionTek has no analytical capability. The VisionTek is also capable of high resolution scanning of the whole slide (WSI) or Partial slide (PSI).

The Sakura VisionTek WSI and PSI capabilities are for research use only in the United States of America and may not be used in diagnostic procedures.

General Description

The VisionTek is an automated bright field digital microscope for histology and cytology. Sample specimens are formalin-fixed, paraffin-embedded (FFPE) tissue on slides, frozen sections, or cytology preparations with glass or film cover slip.

The VisionTek operates interactively, similar to a standard light microscope, with the added benefits of digital technology. It can also operate in automated whole-slide imaging mode.

Upon insertion of the slides, the VisionTek acquires an overview image of each inserted slide, including imaging of their labels. Overview images of the slides are displayed, allowing the user to navigate to positions within the specimen. One can then zoom in to view high resolution images "live" of selected areas. At the highest resolution, sub-cellular and sub-nuclear features can be analyzed.

The instrument has the following features:

- Live examination of specimens
 - Panning scrolling zooming focusing
 - High resolution scans of selected window
- Automated partial or whole-slide scanner
 - Walk-away operation
- Viewer allows examination of stored images
- Slide Tray allows examination and comparison of up to 4 slides simultaneously
- Save snapshots for instant communication
 - single or multiple view
- Annotations can be added to live views and saved images
- Smooth, seamless magnification up to 40x:
 - Auto-switching of three optical magnifications
 - Automatic digital zooming
- Fast and reliable autofocus
- Manual focusing in live mode
- Store multiple focal planes (Z-stacks)
- Tissue recognition assists scanning and focusing

Safety Precautions

WARNINGS, CAUTIONS and NOTES are provided throughout this manual to indicate levels of potential hazards as defined below:

WARNING

IDENTIFIES A POTENTIAL HAZARD IN WHICH FAILURE TO FOLLOW INSTRUCTIONS MAY RESULT IN SERIOUS INJURY TO THE OPERATOR AND/OR OTHER PERSONNEL.

CAUTION

Indicates a potential hazard in which failure to follow instruction may result in damage to the instrument and/or other property, or may result in poor performance.

INTRODUCTION

NOTE

Indicates a reminder or other helpful information.

VisionTek Cautionary Notes

CAUTION: Do not remove or change objectives. The installed objectives are calibrated.

CAUTION: Ensure that there is enough space to the rear of the instrument for proper operation and airflow.

CAUTION: Do not remove the housing of the instrument.

CAUTION: Do not attempt any repairs or alterations.

Physical Characteristics

Cabinet Front (Figure 1-A)

Keypad and LED Indicators

Provide the operating controls and indicators necessary to operate the instrument and monitor system status.

The Keypad buttons provide the following functions:

POWER

Press for >5 seconds to turn on the instrument. (The main switch rear right panel must be ON).

EJECT

Press to eject the Slide Tray. (Only valid when the software is not running).

SCAN

<currently inactive>

The LED Indicators provide the following information:

POWER (Green)

Light is ON only when the power button has been pressed and the instrument is running

BUSY (Orange)

Light ON – Action in process

ERROR (Red)

Light ON – Error state

Slide Tray Door

Provides access for inserting and removing the Slide Tray with specimens.


Figure 1-A

Delivery Contents

VisionTek Digital Microscope includes the following:

- 2 Slide Trays
- 2 x Cable, USB
- Cable, Gigabit Ethernet
- Power cord, 110V, North America NEMA 1-15 or 240V SCHUKO

Additional hardware included:

PC:

- Windows 7 x64 Ultimate SP1 English
- Intel Core i7-2600 (or higher)
- 8 GByte DDR-3
- 120GB primary HD+ 1000 GB secondary HD, 3,5inch,7200rpm

- 6x USB2
- DVI or Display Port (DP) monitor output,
- RS-232 Serial
- 1000 Base-T Ethernet
- Intel Pro/1000 PT Server network card
- QWERTY keyboard and mouse

Monitor:

- TFT monitor with S-IPS or PVA panel
- 24 inch
- 1920 x 1200 pixel
- Display Port or DVI
- Pre-calibrated

Slide Tray

Power Cord

Gig-E camera

2x USB

Figure 1-B


INTRODUCTION

Cabinet Rear (Figure 1-C)

Power Switch — connects the instrument to the main power. Switch ON (I) to use the instrument.

Power Cord Connector — accepts standard power cords for all countries via the IEC-320 plug.

USB, Gig-E Connectors — accepts the provided cables for connection to the PC.

PC Rear (Figure 1-D)

Power Switch — connects the instrument to the main power. Fuses are located in this power entry module.

Power Cord Connector — accepts standard power cords for all countries via the IEC-320 plug.

USB, Gig-E Connectors — accepts the provided cables for connection to the VisionTek.

IMPORTANT NOTES:

- Plug VisionTek USB camera to the second row alone, in the right port.
- Connect VisionTek USB control to the top row, in the far left port.
- The two USB ports connected to the VisionTek need to be diametrically opposed in the PC ports.
- Plug keyboard and mouse into the USB hub of the monitor.
- Connect bottom port to the VisionTek Gig-E Adapter Card.
- This is the built-in network regular LAN connection of the PC. It is NOT to be connected in any way to VisionTek.


Figure 1-D

INTRODUCTION

Specifications

Models Covered By This Manual

Product Code	Name/Description
9000	Vision-Tek® 115VAC – USA
9002	Vision-Tek® 230VAC – Europe
9004	Vision-Tek® 100VAC – Asia

Power Requirements

85 to 264VAC, 47-63Hz, single-phase power cord connector:
IEC/EN 60320-1/C14

NOTE: Universal power supply.

Power Ratings

100V operation: 4 Amps

240V operation: 2 Amps

Environmental Requirements

Operating:

Temperature Range: +15°C to +35°C

Relative Humidity: 10% to 75%, non-condensing

Storage:

Temperature Range: +5°C to +40°C

Relative Humidity: 5% to 85%, non-condensing

Dimensions

Height: 60 cm (23.6 in.)

Width: 50 cm (19.7 in.)

Depth: 50 cm (19.7 in.)

Weight: 35 kg (77.2 lbs)

Instrument Life Expectancy

7 years (with instrument powered on 24 hours a day,
every day)

Safety Standards

Electrical Safety:

For USA:

UL 61010-1 (2nd Edition)

CAN/CSA-C22.2 No. 61010.1-04

For Europe:

IEC 61010-1 (2nd Edition)

IEC 61010-2-101

Regulatory Information:

CE Certified in compliance with the requirements of EMC
Directive, LV Directive in Europe (Europe)

IEC 61326-1:2006

This device is tested to comply with FCC Rules Part 15.

INSTALLATION

General Information

This section provides information on determining a location for, and installing, the VisionTek® instrument. Installation should be performed by a qualified instrument service representative *only*. The VisionTek instrument must be installed correctly to ensure proper operation and service.

Read this operating manual carefully before attempting to operate the VisionTek instrument. Follow all instructions carefully.

CAUTION: The VisionTek is a precision instrument and must be handled accordingly. Rough handling or dropping will disturb or damage internal components. Always handle the instrument with care.

Environmental Factors

As with all sensitive electronic equipment, avoid prolonged exposure to excessive temperature and humidity. Temperature and humidity should be held relatively constant to obtain the highest degree of operating stability. The ambient operating temperature range of the instrument is 15° to 35°C (59° to 95°F). The ambient operating humidity range is 10% to 75% relative humidity, non-condensing.

Locate the instrument in a **well-ventilated** area, avoiding exposure to corrosive vapors, direct air currents, or temperature extremes. Avoid proximity to direct sunlight, open windows, sinks, ovens, open flames, hot plates, radiators, and dry ice baths. Locate the instrument away from any equipment that consumes a high voltage or large current, including large refrigerators and ovens.

Do not cover or block the air intake of the fan. Maintain enough space to other objects behind the instrument (about 4 inches / 100mm), to allow adequate airflow.

Since the VisionTek is a precision optical instrument, it must be installed on a very stable, solid tabletop. The environment should be free from vibration and shock.

Be sure the instrument is located near a power source that meets the electrical requirements (voltage and amperage) specified on the rating label located on the rear of the instrument. The power receptacle must be grounded and should be a clean, noise-free, dedicated line.

Ensure that the selected installation site provides sufficient clearance to allow for proper operation of the instrument, and to provide adequate ventilation.

Unpacking

WARNING: USE CAUTION WHEN LIFTING THE INSTRUMENT. AT LEAST TWO PERSONS SHOULD LIFT, OR WHEN IN DOUBT, USE MECHANICAL ASSISTANCE.

CAUTION: Always lift the instrument carefully by its base when removing from the packaging—see Fig 2-A. Do not put excessive force elsewhere on the housing. Avoid shock to the instrument.

Figure 2-A

Unpack the delivered contents carefully, and retain the packing materials for future transport. Station the instrument on a solid tabletop, with the PC preferably on the floor below and monitor nearby.

WARNING: USE CAUTION WHEN LIFTING THE INSTRUMENT. AT LEAST TWO PERSONS SHOULD LIFT, OR USE MECHANICAL ASSISTANCE.

(See attached document "SVDM01-D1-0003 Uncrating Instructions").

INSTALLATION AND SETUP

Unpacking and Installation

See document "SVDM01-D1-0002 Installation Manual".

Delivery Contents

Description	Quantity
VisionTek, including cables	1
PC	1
Monitor	1
Mouse	1
Keyboard	1
Accessories: Slide Tray	2
Operating Manual	1

Work Surface

The VisionTek is a very sensitive, high resolution microscope, and must be installed on a stable, vibration-free work surface. For example, do not use tables which are on casters.

Warnings and Precautions

WARNING: DO NOT REMOVE THE INSTRUMENT HOUSING. THIS IS ONLY TO BE DONE BY TRAINED SERVICE PERSONNEL.

WARNING: DO NOT REACH INTO THE DOOR OF THE INSTRUMENT WHEN IT IS IN OPERATION.

For optimal results when operating this precision instrument, please follow the suggested operational details below:

- Avoid dust, fingerprints, and contamination on the slides.
- Use slides and materials which are of high quality standards.
- Ensure that slide preparation and staining methods are acceptable.

OPERATING INSTRUCTIONS

Operating Instructions

System Start

The order in which the system is turned on is important. Turn on VisionTek® instrument by switching the ON/OFF switch in the back of the instrument on the lower right side. Then press the [POWER] Button (1) on front panel and hold for two (2) seconds.

- System will perform an internal validation drive: "Busy" light ON – (2).
- System is ready to use when green LED is turned ON, and Busy light is OFF.


Figure 1A


Figure 1B

Turn on system PC, launch Microsoft Windows to launch the VisionTek operating software, and click on the appropriate icon as seen in Figure 2.



Figure 2

OPERATING INSTRUCTIONS



System Operation

Slide Loading

Slides are loaded onto the Slide Tray as shown in Figure 3A. Position the Slide Tray so that the handle is facing user. The slides will then be placed on the tray with the labels up and on the left of the tray. Slides should be properly inserted into their intended slot before loading the Slide Tray into the VisionTek.

The Slide Tray is inserted into the VisionTek using the track as a guide. Once inserted, system operations and pre-scan are activated by clicking the tray into the track (needs a firm push and a “click” is heard).

Slide Tray Loading

The Slide Tray is loaded into the front slot of the system. Press the [EJECT/INSERT] button on the main screen of the software to open the door. The track for the Slide Tray will automatically present for the Slide Tray insertion. If software has not been launched, open the door by pressing on the front panel [EJECT] button.

 The [EJECT] button on the front panel is inactive if VisionTek application is running, to prevent wrongful ejection of the tray during automated scan.

Software Operation for Single Live Slide Analysis

Startup Screen

When the software is launched, the main screen offers three options on the top menu bar left: MICROSCOPE, SCANNER, or SETTINGS. If a Slide Tray is inserted, the instrument will automatically default to Microscope Mode. Just by clicking on the appropriate button, the options of Scanner or Settings can be selected.

Figure 5

Upper Tool Bar

At any time during the process, it is possible to switch between the different modes. The current process will be aborted but there will be a warning to save current data.

To obtain information about the system, place the cursor over the Sakura VisionTek logo in the upper left corner of the screen. The information that appears is shown in Figure 6.


Figure 6

The opposite upper right corner displays the help menu icon. Access to an online help source is located here. Minimize or close it using the standard Windows® "—" or "X".


Figure 7

Startup screen: Microscope Mode

When [MICROSCOPE] Mode is selected either by default or selection, there is an option for:

- Live Slides in the Slide Tray.
- Digital Images only, previously recorded.
- The "Eject/Insert" to remove the Slide Tray.

Figure 8

Slide Selection

When ready to operate the instrument in Microscope mode, the Slide Tray is inserted into the VisionTek with a firm push until click sound heard, then the system operations and pre-scan will automatically be activated.

While VisionTek is pre-scanning the slides, a wait screen will be displayed:


Figure 9

OPERATING INSTRUCTIONS

Next, the instrument will display an overview of all the slides currently on the Slide Tray. These slides are readily available for selection and review as soon as they are visible (Figure 10A). It is possible to allow all the slides to be visible in the overview before selection (Figure 10B).

Figure 10A


Figure 10B

To select a slide, click the box to the left of the slide and then click [OK]. Once the slides have been selected, the instrument provides another view of just the slide(s) selected.


Figure 11A

Figure 11B

Slide Navigation – Description of Mouse Wheel Settings

Navigation within the slide may be done in 4 unique ways:

- Press down on the mouse roller bar and slowly set the pan speed by moving the icon in the direction desired to navigate. This activates the internet scroll function.

Figure 12

- Click and hold the left side of the mouse to grab and move the section in the desired direction. The cursor will turn to the hand icon.

Figure 13

- Move the green navigation indicator in the navigation window located in the upper right corner of the screen.

Figure 14

- Use the keyboard arrow keys to move tile by tile in the direction of choice.

Figure 15

Changing magnification of the view and focusing

To change the magnification for the view, there are two options:

- Use the mouse
 - Left Mouse click: Zoom into image (increase optical lens magnification)
 - Right Mouse click: Zoom out of image (decrease optical lens magnification)
 - Mouse wheel: Adjust image focus live (this feature is not available for overview)

Or

Use the magnification pull down menu

This icon displays the current magnification setting. Magnification changes may be made either by a left/right mouse click or by selecting the appropriate magnification from the pull down menu of this icon.

NOTE: when zooming in to increase magnification, the instrument will automatically center the area that has been selected for the magnification view.


Figure 16

OPERATING INSTRUCTIONS

Full Screen Icon

This icon allows the user to maximize the image display to full screen. When activating to full screen, hide all other views and the current view will fill the screen. The icons on the right toolbar will remain accessible, allowing the user to work with the chosen view in a comfortable manner.

Minimize

This icon minimizes (or hide) the current scene.

NOTE: To retrieve this hidden scene, open the Slide Tray explorer and retrieve the scene from the Slide Tray.

Close Scene

This icon closes the scene. The scene will be permanently removed from the current Slide Tray explorer, unless the slide is already pinned. Pinned Slides will be moved into the Digital Slide category of the tray, and automatically restored in the next session.

Identification Bar

This color coded bar indicates whether the current scene is live or digitized. When the bar is green, the scene is live. Any other color is an indication of a digitized image. The colors change according to the number of digitized scenes available for one slide.

Slide ID and Scene Number

The Slide ID currently displayed is the position of the slide viewed in the Slide Tray (top =1 and bottom =4). The slide ID may be customized by typing directly into the Identification Bar field in the Slide Tray.

Navigation Overview

The Navigation overview is always displayed in the upper right corner of the image. This overview allows the user to physically locate the coordinates of the objective on the microscopic slide. This position is recorded on live views as well as digitized multi-scene views.

Insert/Eject

This icon ejects or inserts a Slide Tray from within the Microscope Live Slides Screen.


NOTE: All live views will be permanently discarded upon Slide Tray ejection since those slides are no longer visible. Digital Slides will be closed and stored in the Gallery. Pinned Digital data will remain available to the user and

will be reopened in the Digital Slide section in next session.

Previous/Next slide position arrow

These two icons allow the user to move to the next or previous slide positioned on the Slide Tray.

NOTE: When the last or slide #4 is reached, by default, the back or next arrow (pointing right) is unavailable. This setting may be changed in the "Quick settings" menu by selecting "Continuous Viewing".

“Add Live Scene” Icon

This is a unique feature of VisionTek. Copies can be made of the live image to allow subsequent comparison with a different area of the same slide. A perfect application for this feature would be viewing different levels of the same tissue on the same slide all at once.

Add Live Scene may be used repeatedly and up to 16 scenes of the slide can be displayed at the same time.

This feature is especially useful when analyzing a tissue microarray, as it allows the comparison of multiple cores, live, on a single screen (with focusing capabilities) – see example below:



Figure 17B

Snapshot Tool

Figure 18

Snapshot Tool

This command allows the user to capture the Field-Of-View (FOV) displayed on the screen (or the active FOV in a multi-scene/slide setting). After activating the command, a classical “Save As” Windows is displayed, allowing the user to define a file name for the recorded picture, and set or accept the default saving path/folder.

There are two options for this icon; single snapshot or multiple snapshots of all the displayed images together. The overlaying icon always defaults back to the last icon/feature used.

Figure 19

OPERATING INSTRUCTIONS

Measurement Tools

When selecting this icon, the user will find a set of tools to use to measure distances or areas currently visible. These tools may be used on a live slide/scene or a captured image. Select the measurement tool to be used by clicking onto it.

Figure 20

Once the tool is selected, drag the cursor to the screen, and draw a line (or square/circle) from the beginning to the end point (this function is similar to the drawing functions classically used in Microsoft Office). Once the object has been inserted (circle), it is automatically deselected. To re-select the object, please click [shift + left click].


Figure 21

With the object selected and the cursor pointed to the object, hold the shift key and right click on the mouse; a context menu will open, allowing the user to specify annotation features. To open the context menu on a previously made annotation, first reselect the annotation by clicking on it ([Shift + Left Click]), or select the specific annotation through the Slide Tray.


Figure 22

Measurement tools can be used in a similar manner for measurement of a round or a rectangular area. Measurement tools may be used on an image in overview as well as on a higher magnification.


Figure 23

OPERATING INSTRUCTIONS

If desired, the measurement performed on a live slide (with the exception of Overview) may be saved digitally, using the single tile scanning function "Mini Scan". Once saved, the corresponding image is displayed in a split window, and new tools become available to further annotate the digital image. The image may be permanently saved using the SAVE function of the Slide Tray explorer.

Figure 24A


Figure 24B


Figure 24C

Quick Settings

The quick settings option allows the user to manage the slide label display (hidden/unhidden/label size), while viewing live slides. There is also a choice to set a continuous slide viewing. This feature allows the user to move back to original position on the Slide Tray after reaching the end of the slide set (go from slide 4 back to slide 1.). See Figure 25 A below.

The Quick Settings icon is located on the bottom of the right menu bar (see Figure 24C)

Image adjustment

The user has the option of adjusting the image "look". The usual image management tools (brightness & contrast), and three different sharpness filters will enhance some image details.


Figure 25A

Figure 25B

Scale

The scale legend as well as the total physical surface measurement as displayed on the screen may be found in the bottom right corner of the image (Figure 26). When the scale legend icon is highlighted, area measurement is displayed by hovering over bottom right end of the image.


Figure 26

OPERATING INSTRUCTIONS

VisionTek Operation for Multi-slide Analysis

Description

Icon Key


Figure 16


Figure 27


Figure 28

"Add Live Slide" Icon

The VisionTek allows to the user to view up to 4 live microscope slides side by side. This feature may be useful when comparing two different stains of consecutive sections of the same case (EX: H&E and IHC).

To do this, click on "Add Live Slide" Icon while viewing the current slide (Figure 27).

The system will bring up an overview of the Slide Tray loaded into the instrument (Figure 28).

"Slide is Already Selected" corresponds to the slide currently being viewed.

Then select one (or multiple) of the other slides to view concomitantly.

Click on the corresponding checkbox, and then select "OK".

VisionTek will open a new window showing an overview of new slide selected.

Navigate on the new slide as described previously, as well as navigate back and forth between all displayed live slides.


Figure 29A


Figure 29B

Concomitant use of “Add Live Scene” and “Add Live Slide”

NOTE: These two features can be used concomitantly.

When comparing two slides, multiple scenes of the first slide (EX: H&E) may be compared with multiple scenes of the second slide (EX: IHC). In Figure 30 below, there are two staining methods and different scenes of the same slides (different sections, note the navigation overview).


Figure 30

OPERATING INSTRUCTIONS

"Add Digital Slide" Icon

Displays views of Live Slides (or Scenes) as well as previously recorded digitized images side by side.

In order to display a digital image concomitantly to a Live Slide/Scene, click on the corresponding icon (see Figure 31).

Figure 31

VisionTek will open the Digital slide gallery for selection. The slide image(s) to be opened are selected and the select box checked (1). See Figure 32.

By clicking on the + sign on the first column (2), a preview of the slide will be displayed. Once all digital slides needed have been selected, click open (3). See Figure 32.


Figure 32

The selected Digital Image will be displayed concomitantly to the Live Slide(s) Scene(s).

NOTE: The live scene in Figure 33 is the top slide (green bar) and the digitized slide is the bottom slide with mauve bar.


Figure 33

NOTE: Important Information: all "Multi-Scene" features may be used concomitantly. It is possible to display:

- Multiple Live slides together with
- Multiple Scenes from a Single Live Slide
- Multiple Digital Images with
- Multiple Digital Scenes from a Single Digital Slide

Currently, VisionTek allows the concomitant display of 16 images mixed of either or live and digitized.

"Add Digital Scene" icon

This feature is equivalent to the "Add Live Scene", but applies to adding an additional digitized image from the active digitized image.

While reviewing a slide image, a scene can be "frozen" to allow subsequent comparison with a different area of the same digital slide. The details are explained below.

While reviewing a slide in the Live Mode, click on the "Add Digital Scene" icon (Figure 34).

Figure 34

A second image will display on the screen.

Navigate and zoom in this second "scene" as described previously.

This feature will allow to the user to compare two or multiple areas of the same scan. A good example of an application for this feature would be comparing cores on the same tissue array.


Figure 35

Digitized Image Files

Saving images as .jpeg or .bmp files

There is the option of saving the live image displayed on the screen as a standard image file for future use in a different application (i.e. PowerPoint). Two unique methods are available:

- Single snapshot
- Multi-View snapshot

Single Snapshot

This method gives a simple snapshot of the current live image scene regardless of what other images might be open on your screen.


Figure 36

OPERATING INSTRUCTIONS

Digital Recordings of Images

NOTE: The snapshot tool: saves images as .jpeg, .tiff, .png, or .bmp files for use in different applications.

The user saves the live scene displayed on the screen as a regular image file, under a Windows Explorer folder, for future use in a different application (i.e. PowerPoint). The file created is of standard size and resolution (order of magnitude: one megabyte)

To do so, user may use snapshot tool.

Figure 37A

This icon will trigger the opening of a dialog box, allowing the user to define the details to be captured in the standard resolution image to be created. An image preview in the dialog box will display the image file as it is being captured.

The user may capture the single active view, or the multi-view of all images displayed on the monitor.

Single View Capture

In this configuration (1), the active view (2) is captured and the preview of the file is displayed in the dialog box (3).


Figure 37B

Multi-View Capture

The snapshot mode can be changed by toggling on the corresponding icon (clicking anywhere in the area will toggle the mode)


Figure 38A

The multi-view snapshot mode will capture all slide views displayed on the monitor.


Figure 38B

OPERATING INSTRUCTIONS

High Resolution Scan

In addition to the snapshots explained above VisionTek is capable of high resolution scans. There are two methods available for scanning the whole slide, partial slide, or selective area on the slide.

- Scan as needed while viewing
- Automated scan for review later

While reviewing a slide image in Live session, a high resolution scan may be desired. When selecting "scan" the parameters for the scan are set by the user.. Select the scan icon from the right menu bar (Figure 39).


Figure 39

Once scan is selected the image will revert back to the overview, so that you can define the scan area. The red horizontal and vertical lines intersect on the last scene you viewed. Simply grab the square and adjust to cover the area to scan. This feature allows the user to scan only the area of interest instead of a predefined size, shape, or area. The user controls it. Once the area has been set, press [OK] (yellow arrow in Figure 40).


Figure 40

Next, VisionTek will return the view of the live scene, at 20x magnification. The user is asked to adjust the focus to ensure a crisp image. Use the roller ball on the mouse to adjust. Once the focal plane is properly set, press [OK] (yellow arrow in Figure 41).


Figure 41

A last menu section unfolds, allowing the user to define:

- The scan resolution
- The number of Z-stack planes
- The distance to be set between the planes

Figure 42

Once all scan parameters have been set, press [START SCAN].

The VisionTek will display the digital image as it is built. During the scan process:

- Zoom in at any time to verify the quality of the image (Figure 43 1);
- Enter text in the comment section and define a unique name for the scan(2),
- Or the scan may be cancelled (3).


Figure 43

OPERATING INSTRUCTIONS

Once the scan is completed, the recorded image will be displayed in a new window (Figure 44), allowing the operator to add further annotations (measurements, text, etc) on the image.

To be permanently saved, the scan and the corresponding slide information must be saved in the Slide Tray explorer (2).


Figure 44


Figure 45

High Resolution Mini Scan

If only field of views (FOV) are necessary to document specifics on a slide, the VisionTek allows the user to perform a High Resolution Mini Scan, saving only the area or FOV displayed on the screen. This is an image equivalent to a snapshot, but with high resolution (Figure 45). This image would be added to the slide file folder, and stored in the Digital Gallery.

When activating this feature, a dialog window will pop-up. Focus and Z-stack planes and z- step size are pre-set by the view on the screen. The user only needs to enter a scan name and a scan comment (Figure 46).


Figure 46

Automated Scanning

When using the automated scanning function, the live image is bypassed.

To proceed with an automated scanning, insert one to four slides (selected for scanning) into the Slide Tray, and insert the Slide Tray into the instrument.

On the startup screen select the Scanner mode (1). See Figure 47.





On the subsequent screen, select the scan profile required for the loaded slides (2) See Figure 48. Multiple slide profiles may be defined in the Settings.

Start the scanning process by selecting [SCAN SLIDES] (3). See Figure 48.

VisionTek will proceed with detecting the presence and position of all the slides.

During the scan process, the status is displayed on the screen:

- The system will first proceed with a pre-scan of the label and the slide (Figure 49A), and define a focus map with various focus points (arrows), the focus points being distributed randomly over detected tissue present in the scan area (Figure 49B).
- The system will then perform an autofocus check for each point of the focal map:
 - A successful focus point is displayed with a green dot.
 - A failed focus point is displayed with a red dot.

OPERATING INSTRUCTIONS

Viewing of Digital Slide

Opening of Digital Slides

View previously recorded digitized images, without any Slide Tray loading requirement. In this case, the software will act as classic "Image Viewer".

Figure 50

Selection of Digitized Images

After entering the digital slide mode, the system will display the Digital Slide Gallery, allowing for selection of one or multiple images to be displayed on the screen. Select by clicking on the select box and then press [Open].

Figure 51

Features


Figure 52A

Same Features as Live Mode

The features available in the Digital slide mode are identical to the features available in the Live Slide Mode, with the exception of the "Next/Previous" feature which moves from one physical slide to the next. Observe in Figure 52A above, those icons are missing (circle).

Additional Annotations

All Digital Slides may be annotated using the measurement tools, and yet there are additional drawing tools (bracket in Figure 52A).

NOTE: When the text tool is selected, text may be entered directly into a dedicated text field. This is done by opening the content menu (Shift + Left Click to select, Shift + Right Click to open content menu).

Addition of Live Slides present in the tray

Live slides can be added to the displayed scene, by using the regular "Add Live Slide" icon. The instrument will then continue to perform the standard preview of labels and slides. Proceed with standard slide selection.

Scan information

As for a Live Slide, the Digital Image displays a scale on the bottom right end of the image, and the corresponding area surface.

Additionally, the digital slide displays the scan resolution used to record the image on the bottom left corner of the slide. In both areas, the information will be highlighted when the cursor hovers over that area.


Figure 52B

OPERATING INSTRUCTIONS

Slide Tray – Slide and Data Management


Figure 53

The Slide Tray explorer is an important piece of the VisionTek operating system, allowing the user to easily manage all the views, annotations, or any information related to the slides currently viewed. It serves the same purpose as a desktop slide folder and will give the operator a good overview of all slides needed to diagnose the particular case.

The Slide Tray is folded in by default. To have the Slide Tray unfold, simply hover the cursor over the right side tool bar at the Slide Tray tab. The Slide Tray will automatically fold open as you move the cursor away. It is also possible to pin the Slide Tray in place so that it will stay unfolded. Use the Pin icon in the upper left corner of the Slide Tray window (Figure 54 circle).



Figure 54

General Slide Tray Features (Figure 55)

OPERATING INSTRUCTIONS

Feature Descriptions (Figure 55)

1. Microscope Live Slide Section

This section of the Slide Tray will list all Live slides physically present in the Slide Tray for managing all corresponding live views, and/or digital images and annotations freshly saved during the Microscope mode.

2. Microscope Digital Slide Section

This section of the Slide Tray will list all Digital slides open, for managing all corresponding digital scenes, and their corresponding annotations and/or comments.

3. Slide ID

This is the default according to the slide position in the Slide Tray (top=1 and bottom = 4), this field may be edited just by clicking directly into the field. This feature allows the user to uniquely name the slide, and therefore facilitate its identification in the digital gallery.

4. Slide Label Image and slide overview

Each section of the Slide Tray is identified by the display of the slide label, and the corresponding slide overview image allowing the user to properly identify the slide and all the current views. Furthermore, this area of the Slide Tray serves two additional purposes:

- Orientation of the slide label image can be temporarily changed during viewing from the default settings, by clicking on the label in the open Slide Tray: each mouse click will rotate the label image by 45° to the right in the Slide Tray as well as in the main scene.
- All views related to this slide may be displayed or hidden (fold/unfold function) by simply clicking on the image overview


Figure 56


Figure 57

5. Slide Scenes

All scenes related to a single slide are listed and displayed under the corresponding slide label and overview. Each scene is underlined with a tool bar displaying the scene number.

There is a hide icon to hide the corresponding scene in the main view; when the image is hidden, this icon will automatically toggle to “unhide” icon.

Figure 58

Click on a specific thumbnail and this scene will become active on the main screen and will be enlarged in the Slide Tray fold, to reflect the working scene. The active scene on the screen is then outlined by a blue background (Figure 59, circles).

To make annotations at this point, use the “full screen” icon in the desired scene to enlarge the image to edit.


Figure 59

OPERATING INSTRUCTIONS

The scenes are color coded:

- All scenes of a live slide are displayed with a green shade background
- All scenes of a digital image are color coded using any other color.
- Each new scan is coded with a unique color, and each additional scene for this same scan will be coded with the same color.

In the example below, this slide has 9 different scans linked to it.

There are three scenes opened for the light pink scan, and two scenes open for the purple scan.


Figure 60


Figure 61

Scenes will be grouped under the slide overhead either/or

- in a checkerboard display if the slide is inactive
- in a linear way for an active slide; to scroll through an active slide scenes, please use the mouse wheel

11. Fold /unfold all scene views

To facilitate the management of the Slide Tray folder, these icons will help to expand or collapse all views present in the Slide Tray, independently of their individual status.

6. Close Slide

This icon closes a slide with all its corresponding scenes with one mouse click. Should the slide be a digital slide and some modifications left unsaved, a warning message will trigger the user to save remaining changes, if desired.

7. Save slide

This icon saves all changes made for a slide in a single mouse click (saves newly created scans, annotations...).

8. Hide/Unhide all

This icon hides or unhides all scenes of a given slide, with a single mouse click.

9. Delete Slide

This icon will delete a given digital slide, with all its related scans, in a single mouse click. All data will be permanently removed from the database.

NOTE: To remove a single scan within a digital slide, user may trigger the corresponding specific icon.

Figure 62

10. Diagnosis Code

Optional field to enter the diagnosis code for the case, once diagnosis is finalized. Diagnosis code is subsequently displayed in the database and allows easy filtering and case search.

12. Virtual Case Session

The virtual case session allows the reviewing of a complex case of multiple slides, and records all interesting fields of views for final review. The specific functioning of this mode is described in a separate paragraph.

13. Close All

The [CLOSE ALL] button saves all digital images with open (unsaved) modifications all at once, (rather than saving all slides individually) before closing the Slide Tray folder and returning to the microscope mode.

OPERATING INSTRUCTIONS

Particular details of a single Microscope Slide

Annotation and Measurement

In each slide scene (whether live or digital), if an annotation and/or a measurement are added to a particular image, the inserted object is listed in the legend field for that image.


Figure 63A

A single click on the object in the legend field easily selects the object in the image (i.e. to modify the properties in the content menu) (hold shift + right click).

A double click on the object in the legend field re-centers the object in the display (i.e. confirming the position).

Hovering the cursor over the object in the legend field will display the name set for the object in the content menu (see Figure 63B).


Figure 63B

Hiding/Unhiding Scenes


Figure 64A

If used, this icon automatically toggles to its "Show" counterpart.

The Hide icon hides a scene on screen, from the Slide Tray.


Figure 64B

The Close icon linked to each slide removes the corresponding slide from the Slide Tray.

If needed, the slide may be added again by simply using the "Add Live Slide" or "Add Digital Slide".

The Close icon linked to each scene permanently closes the corresponding scene in the main view, but keeps the slide active in the Slide Tray.

OPERATING INSTRUCTIONS

Digital data recorded during Live Microscope Mode

When digital data is recorded during the Microscope Live mode, a [SAVE] button will appear in the Slide Tray toolbox of the corresponding slide. If clicked, a dialog box appears for confirmation or re-edit.



Figure 65A

Particular details of Digital Slide Section

The Slide Tray waste basket displayed allows the user to delete an entire slide from the database, independently of the number of scans attached to it.

The waste basket displayed in the Scan legend area allows the user to delete this specific scan, but keep all other meta information and additional scans available in the Digital Gallery.

Figure 65B


Figure 66A


Figure 66B

By clicking on the scan scene in the Slide Tray, the system will toggle to text field area and display comments entered during Live scan, or a default comment associated during auto scan process. The same text field area will display some detailed scan information such as scan date, area, number of Z- stacks, etc.

To revert back to scene thumbnail, simply click on the text area again.

The Virtual Case Session

When reviewing multiple slides from a complex case on a manual microscope, a pathologist has to memorize the areas of interest on the various slides or mark areas using a marker, so that he/she may revisit these areas of interest, to confirm his diagnosis hypothesis.

The Virtual Case Session enables the user to go through a similar process and capture all areas of interest, during the case review. The collection of snapshots so created may then be reviewed upon case completion.

Valuable slides showing areas of diagnostic interest may then further be saved in the databank, and be assigned the proper diagnosis code, for further retrieval.

The virtual case session may be activated by clicking on the appropriate icon. A virtual case session in progress is symbolized by a green light activated in the given icon.


Figure 67

OPERATING INSTRUCTIONS

By starting a virtual case session, the user will slightly change the workflow of the data capture.

- While reviewing the first slide of the case, capture all areas of interest using the "Mini Scan" functionality. A preview of the captured area will appear as usually, but once the OK button is clicked, the snapshot is automatically saved and minimized, allowing the user to seamlessly navigate through the slides.


Figure 68

The user may review all first four slides of the case in such a way.

In the example below, the 4th slide of the case is currently being viewed, and several scenes of slide 1 through 3 of this case have been captured and automatically saved.

Once the first set of four slides of the case have been reviewed, proceed with ejecting the tray and loading the four subsequent slides. Capture all interesting scenes, as desirable.


Figure 69

Upon reviewing the last slide, and tray ejection, all saved slides are still available on the Slide Tray. All captured scenes can be unhidden with a single click, for global reviewing.


Figure 70


Figure 71

OPERATING INSTRUCTIONS

Once diagnosis is completed, the user may select to save some or all slides of interest (1), and assign them the proper diagnosis code (2). Slides may also be automatically closed (3).


Figure 72

Digital Slide Gallery				
Open Mode		Export Mode	Delete Mode	Import
	Case ID	Slide ID	Diagnosis Code	Date
☐	1234-8		6300 3526	2/13/2013
☐	1234-5		6300 3526	2/13/2013
☐	1234-4		6300 3526	2/13/2013
☐	1234-2		6300 3526	2/13/2013
☒	1234-1		6300 3526	2/13/2013
☐	Skin Sample		0	2/13/2013

Pin Function


Figure 73


Without pinning the images, ejecting a Slide Tray automatically closes all open Views – Live and Digital – displayed on the main screen, and all corresponding slides are permanently closed.

In some cases, some digital data recorded from a Live Slide as well as Digital slides present on the main screen, may need to be maintained in the main screen, despite ejection of the Slide Tray.

To maintain the images, simply pin the image to the Slide Tray using the Pin icon. All digital views will be closed upon tray ejection, and reopened in the Digital Slide section of the Slide Tray in the following session.

NOTE: Similar to taping a note to a board, a pinned slide will always populate the Digital area of the Slide Tray until the slide is purposely removed from the Slide Tray by using the “Unpin” and/or “Unpin all” function.

The Digital Gallery allows the user to easily manage all data stored in the data bank. Several modes are available, and are selectable through the Tabs at top of the window.

This modus allows the user to open and review one or multiple slides in the data bank. A preview of the slide and its contents may be displayed by unfolding the various information layers of the data bank.

- Number of scans: number of sub-regions digitally recorded for a given slide
- File size: total file size of the digital slide (all sub-regions included)
- User: author of the digital recording (linked to Windows log in).

The first layer in the databank displays all general information related to the digital slide:

- Case ID: currently, it records the slide position in the tray or the ID manually entered by the user. Once barcode recognition is implemented, this field will be automatically populated with the appropriate entry
- Slide ID: once barcode recognition is implemented, will be populated with appropriate entry – should the lab actually distinguish between the various slides of a case
- Diagnosis code: manually entered by pathologist (SNOMED* code for instance)
- Date: date of digital recording
- Scan time: time of digital recording
- SNOMED = Systemized Nomenclature of Medicine

OPERATING INSTRUCTIONS

Second layer

The second layer will display (opens when hovering the cursor over the first layer entry for a particular slide):

- a thumbnail of the entire slide
- a thumbnail of scan region and its individual file size
- scan name and comments entered by the user


Figure 76

Third layer

The third layer will display all annotations and measurements added by the user (if any), and their corresponding values.



Figure 77

Filtering tools

The Digital Gallery is outfitted with various filtering options.

- Displaying only slides with associated scans. The default setting is to display only slides with associated scans. Should this parameter be unchecked, the system will display all slides ever loaded by the instrument, and the preview of those slides may be found in the first layer of the gallery.

Figure 78

- Filtering by column head

The columns may be filtered in ascending or descending order by first selecting the column and then selecting the filtering order.

Figure 79

This filtering process may be applied to the following columns:

- Case ID
- Slide ID
- Diagnosis Code
- Date
- Number of scans
- User

- Advanced filtering options


Figure 80


Figure 81

OPERATING INSTRUCTIONS

The Digital Gallery includes an Advanced Filtering option allowing the user to set multiple filtering parameters (additional or mutually exclusive). The filtering algorithm follows the classical model:


Figure 82

Export/Import mode

Export slides and images scanned on the VisionTek to another user database, or vice versa, import data produced by a different instrument.

Once the Digital Gallery is displayed, select the mode to work with by clicking on the proper tab at top of the window.


Figure 83

Checkmark (1) the slide(s) to be exported

- All scanned sub-regions may be exported at once.
- or
- Some specific sub-regions may be exported individually.

Example 1: the entire slide, with its 2 sub-regions, will be exported.


Figure 84

Example 2: only the 2nd sub-region will be exported, together with the slide metafile (whole slide image...).


Figure 85

Click on Export at bottom of screen (circle)


Figure 86

OPERATING INSTRUCTIONS

- A window will pop up, then select the proper folder (or external drive) to export to.


Figure 87

- Hit OK.
- A window will display the export process progress:
 - The Overall bar summarizes the overall export time to expect, should multiple slides have been selected.
 - The Current bar displays the progress in exporting a specific slide.


Figure 88

IMPORTANT NOTE:

The exported file is identified with a unique slide identifier (GUID).

Should a slide with an identical Case ID be imported into the same database, it would be inserted as a new object (due to a different GUID), and would therefore not overwrite the existing slide identified by an identical Case ID.

Import function

The import function is similar:

- Go to Digital Slide session.
- Click on Import.
- Select the drive or folder to import from.
- Hit OK.
- The imported slide will be inserted into the local database, and displayed on the Digital Gallery list on the next access.


Figure 89

Settings

Predefine some defaults options for tools or features of VisionTek.

General Tab

The General tab allows the user to specify:

- The language to be used.
- The preferred default path to export scanned slides or picture snapshots.
- The preferred default path to export Snapshot images.


Figure 90

Display tab

Display tab - Annotations

The first category of this tab allows the user to define the standard colors, fonts and drawing characteristics.

Display tab - Display

The second category of this tab allows the user to define:

- Slide label visibility.
- The size of the label to be visible by default.

The display mode:

- "Pathology mode" will display the entire area captured by the sensor in a compressed manner in order to display the largest Field of View possible. This allows a fast slide analysis without any loss of quality.
- "1:1" will display each pixel from the sensor to one pixel of the screen. Therefore, the amount of slide area displayed will be restricted (smaller Field of View). In exchange, the user will be given the impression of a higher magnification. This tool may be useful for some applications needing high nuclear resolution

OPERATING INSTRUCTIONS

The difference between the two modes becomes obvious if the user reviews some screen shots of the same object for the two modes.


Figure 91

NOTE: The 1:1 mode favors the resolution and always displays 1 pixel to 1 pixel. Therefore, in a split window, the Field of View will be cut in half for the 1:1 mode, while the image will be compressed to preserve its area integrity in the pathology mode.

The cursor key overlap defines the amount overlapping Field of View pixels when navigating through the tissue using the arrows on the keyboard (screening mode). Scroll sensitivity defines the speed of image displacement if activating the scrolling mode.

Default Label Rotation angle defines the default angle to be used in the Scene screen (Vertical or horizontal). By clicking on the arrow, a drop down list is displayed allowing to the user to select the default angle to be used.

NOTE: This display angle can also be changed temporarily during a microscope sessions (Live or Digital) by clicking on the label itself in the Slide Tray.

NOTE: These settings can also be adjusted during image viewing by using the Live Image adjustment function.

Figure 92

Display tab - Image Adjustments

The third category of this tab allows the user to define:

- The default brightness and contrast of the displayed images.
- Applying a sharpening filter by default.
- The settings to be applied to the sharpening filter by default.

Figure 93

OPERATING INSTRUCTIONS

Microscope Tab


Figure 94

This tab allows the user to define:

- Continuous slide viewing (i.e.: go back to slide 1 once last slide has been read).
- The starting objective resolution (able to skip the Overview image if deemed necessary).

Scanner Tab

Auto-Scanner Mode

The first category of this tab allows the user to define the automated scanning parameters, by defining scan profiles:

- Existing scan profiles may be selected using the drop down box.
- Existing scan profiles may be deleted using the [REMOVE] button. New scan profiles may be created using the [ADD] button.

General Fields

The scan profile name allows the user to edit the name of the current profile considered.

Two editable fields allow the user to enter a default scan name and short description.


Figure 95

Scan area field

The scan area field allows the user to select between two types of automated scanning mode:

- Automatic tissue detection: The system will locate the tissue on the slide, and scan only areas displaying tissue. Surrounding empty glass area will not be scanned. This feature will provide scan files of minimal size and help reduce the total amount of storage space needed to preserve the data.
- Fixed scanning area: The system will scan a defined area on the slide, whether this area contains visible tissue or not. Visible tissue will be used to determine the autofocus plane as a primary step; should no tissue be visible on the slide, the system will revert to a phase interface autofocus method, as a fallback position. The fixed scanning area is useful to scan cytology spots (density of cells may be low and hardly recognized by the system and they are contained to a specific area).

Automatic Tissue Detection

Working area: defines the area on the slide where VisionTek has to search for tissue. Reducing the area to the acceptable minimum will help reduce “Spam” scan, which may occur if some printing or writing is present on the slide (“SuperFrost +” print, for example). The working area can be set by simply modifying the size of the blue outline (click and drag on rectangle corners). Minimum Object size: defines the minimal possible size of the objects to be detected. Current lower threshold is defined as default: 1x1mm


Figure 96

Fixed Scanning Area

By enlarging or reducing the yellow outline on the slide, the user defines the area to be scanned for all slides (when using that profile), independently of the presence of tissue.

This feature is useful to scan, for instance cytology spreads or tissue arrays. Using the fixed scan area will ensure that the entire TA is properly scanned.


Figure 97

OPERATING INSTRUCTIONS

Auto Scan Settings: Scan resolution and Z-planes

The Resolution drop-down box allows the user to define what default scan resolution should be used for this particular profile; default scan resolution may be different for different applications and needs.

The number of Z-stack planes and distance between them may be set for each profile.

It is possible to define whether the Slide Tray is to be ejected, upon scanning completion.


Figure 98

Live Scan

The second category of the Scan settings tab allows the user to define the default resolution to be used when performing a scan during a live microscope slide session.

Figure 99

Import/Export Function


Export slides and images scanned on the VisionTek to another user database, or vice versa, to import data produced by a different instrument.

Export Data

In order to proceed with an import or export procedure, proceed to the Microscope Mode/Digital Slide menu.

Figure 100

- A window will pop up, then select the proper folder (or external drive) to export the file to.
- Press [OK]. A waiting screen will appear during the export process, and a toolbar at the bottom of the screen displays a progress bar.

NOTE: The exported file is identified with a unique slide identifier. Should this slide mistakenly be imported again into the same unit (database), an error message will be displayed warning that this slide already exists in the database.

The Digital Slide Gallery will be displayed:

- Select (1) the slide to be exported (best to export slides one by one, given the size of the files).
- Click on Export (circle).


Figure 101

OPERATING INSTRUCTIONS

Import Function

The import function is similar:

- Go to Digital Slide session.
- Press [IMPORT].
- Select the drive or folder to import from.
- Press [OK].
- The imported slide will be inserted into the local database, and displayed on the Digital Gallery list the next time it is accessed.

Alt-F4:

- Leave application.

Keyboard Shortcuts

Currently, keyboard shortcuts are set by default in the following manner:

Escape:

- Deselect all annotations in current scan/slide.
- Leave full screen.
- Leave edit mode in slide name control.

Delete:

- Deletes selected annotations in current scan/slide.

Page Up/Down

- Change focus in current scene.

Home/End:

- Zoom in/out in current scene.

Left/Right/Up/Down:

- Move around in current scene.

Shift-Clicks:

- Annotation selection and context menu.

Ctrl-Click:

- Toggle to full screen for current scene.

RECOVERY PROCESS

Overview

This section describes the recovery process for any operation that is suspended due to an error condition. See "Section 8, Troubleshooting" for specific troubleshooting procedures for a given error.

It is important to note that not all conditions causing an error are immediately recoverable. Some error conditions require a service representative to diagnose the problem and may require service attention. In this case, please contact the Sakura Technical Support Department as 800-725-8723, option 2 (U.S. customers only).

Actions for Recovering from an Error Condition

Error message displayed

When the instrument has detected an error during operation, a message window displays automatically. When it appears, *write the message down*, for further communication with Service.

Instrument not powered on

If the VisionTek® application is launched while the instrument is off, an error message will appear and only the View Mode will be available. Turn on the instrument, and re-launch the VisionTek application.

No Slide Tray Inserted

An error message is seen when "Insert Slide Tray" is pressed, and no Slide Tray is inserted. This error may also be seen when the tray is not properly inserted. Listen for the "click".

Insert a Slide Tray with slides to be analyzed.

Software freezes, does not respond to input

a)

- Close the VisionTek Application using the Windows Task Manager.
- Launch the VisionTek Application anew.
- If the condition persists:

b)

- Shut down and re-start Windows.
- Launch the VisionTek Application anew.
- If the condition persists:

c)

- Cycle the power of the VisionTek off then on via the POWER button.
- Launch the VisionTek Application anew.
- If the condition persists:

d)

- Check the connections between the VisionTek and the PC (reference the Troubleshooting section below)

Error message during Automated Scan or HR Scan

(WSI/PSI)

Errors may be caused by contamination on the slide, or if the slide composition/quality is not adequate.

Clean the slides and repeat.

RECOVERY PROCESS

CARE OF THE INSTRUMENT

Maintenance as Needed

Slide Tray Cleaning

Slide Trays must be kept clean to minimize ingress of dust and contamination into the VisionTek®.

Clean by wiping them using a non-abrasive lint-free cloth, if necessary, use a mild detergent. Do not use acids, caustics or chlorinated solvents. Ensure that the seating surface for the slides is free from contamination.

Annual Maintenance

A yearly maintenance by Sakura Technical service personnel is recommended.

Daily Maintenance

Exterior Cleaning

Keep the exterior of the cabinet free of dust and clutter. Generally, wiping with a clean cloth moistened with water is all that is needed to keep the cabinet dust free. Do not use solvents of any kind on the front control panel. If the instrument must be disinfected, wipe the cabinet down with 70% ~ 80% ethanol or isopropanol, and allow to air dry.

CARE OF THE INSTRUMENT

TROUBLESHOOTING

General Troubleshooting Procedures

The following Troubleshooting Guide contains information regarding general issues that could occur during operation of the VisionTek®. Recommended remedies are included.

Operator involvement may be needed in the following two cases:

After connection to a PC, VisionTek does not operate or does not operate correctly.

- Check the connection of all the cables between the VisionTek, the PC, and the monitor.

Service Information

When a Problem Occurs

If you have consulted and followed the instructions in the troubleshooting guide and the issue has not been corrected, call Sakura Technical Support Department for assistance.

Before calling for instrument service, please have **model number, installation date, and warranty period** ready for our Technical Support Representative. For convenience and reference, record this information in the blanks below.

Model Number: _____

Serial Number: _____

Installation Date: _____

Warranty Period: _____

Where to Call for Service

If located within the United States, contact the Technical Support Department of Sakura Finetek USA, Inc. by calling toll free:

1-800-725-8723

In countries other than the United States, contact the nearest authorized VisionTek® instrument distributor or representative for service information and assistance.

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All Rights Reserved