Approved by Director of
TITAN BIOTECH LTD.

SURESH CHAND SINGLA

Instructions for use of the medical product
«NUTRIENT AGAR»,
Produced by TITAN BIOTECH LTD., India

With attached certificates

11 MAR 2017

ATTESTED
NOTARY PUBLIC
DELHI (INDIA)

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№	Parameter	Specification
1	Appearance	
1.1	Physical state	Powder
1.2	Color	from straw-yellow to dark-yellow
1.3	Hygroscopicity	Hygroscopic
2	Dissolubility	28 g/l
3	Transparency and color of the finished medium	
3.1	Transparency of the finished medium	Transparent, with mild opalescence
3.2	Color of the finished medium	Yellow-cream
4	pH of the finished medium at +25°C	7.4 ± 0.2
5	The strength of the jelly of the finished medium	Not less than 800 g/sm ²
6	Gel temperature	30-37°C
7	Melting temperature of the jelly of the finished medium	Not less than 80°C
8	The melting time of the jelly of the finished medium	Not more than 1 hour

Quality control
Sterility testing of the finished nutrient medium
A macroscopic testing should not reveal any microbial growth after incubation at 20-24 ° C and 30-34 ° C for 5 days.

4. ANALYTICAL CHARACTERISTICS

№	Parameter	Specification
1	Determination of specific activity	
1.1	Sensitivity of the medium	The medium should ensure the growth of the test

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Results after aerobic incubation at 35-37 ° C for 24-48 hours

5. PRECAUTIONS

Nutrient agar is intended for in vitro diagnostic use only.

Only for professional use.

Avoid inhalation, contact with eyes, skin, mucous membranes and clothing.

After contact with sensitive areas, immediately rinse with water, wash skin with soap and water.

Avoid microbial contamination of the reagent.

It is strongly recommended that the control strains of microorganisms be potentially infected and observe all necessary precautions when working with them.

Reusable devices must be sterilized in accordance with the accepted procedure after use, but the preferred method is sterilization in an autoclave for 15 minutes at 121 ° C. Disposable items should be sterilized in an autoclave or incinerated. Spilled potentially contaminated materials should be immediately removed with absorbent paper, and contaminated areas should be treated with a standard antibacterial disinfectant.

Do NOT use sodium hypochlorite. Materials used to clean spills, including gloves, should be disposed of as biologically hazardous waste. Do not take the substance into the pipette with your mouth. When handling specimens and analyzing, use disposable gloves and eye protection. At the end of the work, wash your hands thoroughly. Do not use control strains of microorganisms with expired shelf life.

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6. SAFETY REQUIREMENTS

Work in the laboratory should be conducted in compliance with internal standards and requirements in force in the country. When working with a medical device, follow the requirements of ISO 15190: 2003.

7. EQUIPMENT AND MATERIALS

For the analysis, the following types of equipment, measuring instruments and materials are additionally required: Laminar box (biosafety box) class II / type A2 /, inoculation loops / needles for sowing, gas burner for creating aseptic conditions during seeding, sterile test tubes, sterile pipettes, Drygalsky sterile spatula, sterile microbiological tampon, sterile saline, disposable gloves Eye protection during sowing, McFarland 0.5 microfossure suspending standard, incubator incubator, timer, disinfectant, sterilization autoclave, autoclave or dry oven for sterilization of dishes.

The need for using a specific type of equipment, measuring instruments and materials required for work, as well as their performance characteristics, is determined at its discretion by a microbiologist engaged in the formulation of the analysis.

8. ANALYZED SAMPLES

Clinical analytes are samples from the human body, presumably containing microorganisms that do not have specific growth requirements (secretion of outer ear, urine). Agar nutrient can also be used to keep or to check the purity of subcultures of colonies that are undemanding to the composition of the nutrient medium of bacteria grown on other nutrient media from different samples.

9. SELECTION AND TRANSPORT OF ANALYZED SAMPLES

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The sampling and transportation of the test samples is carried out in accordance with the requirements of the Order of the Ministry of Health of the USSR of April 22, 1985 No. 535 and MU 4.2.2039-05. All manipulations with clinical specimens should be carried out in accordance with SP 1.3.2322-08, in order to avoid contamination with microorganisms that may be contained in the samples being tested. It should be remembered that microorganisms can be unevenly contained in the test material, therefore it is very important to take the biomaterial as far as possible from different places of the infection site and in sufficient quantity for analysis.

9.1 Secretion of outer ear

Biomaterial Collection

In infectious diseases of the outer ear, the skin is treated with 70% alcohol, followed by washing with sterile saline. The secretion is collected on a sterile disposable tube probe or tube with a transport medium. The tampon is placed in a test tube and delivered to a laboratory.

Transportation of biomaterial

The obtained samples are immediately delivered to the laboratory, where they are immediately subjected to analysis. If it is not possible to make immediate seeding, it is recommended to store the samples in the refrigerator at a temperature of 2-8 ° C, but not more than 3-4 hours. In the case of using a test tube with a transport medium, the test samples are stored at room temperature (18-20 ° C).

9.2 Urine

Biomaterial Collection

The morning urine is used for microbiological analysis (collection of daily urine for these purposes is unacceptable).

Collection of urine should be done before the start of drug treatment. If the purpose of the study is to assess the effect of the therapy, then urine culture is performed at the end of the course of treatment (on the 5-7th day after the abolition of antibiotics or uroseptics), with the obligatory indication of the drug that was treated. The collection of morning urine is carried out after a careful toilet of the external genitalia and the

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anus area.

In a disposable sterile container with a screw cap,
an average portion of freely released urine is
collected in an amount of 3-7 ml.

Transportation of biomaterial

Delivery of material to the laboratory is carried
out within 1.5-2 hours after collection. It is
allowed to store urine in the refrigerator (at a
temperature of 2-4 ° C), but not more than 3-4
hours. In the case of delivery of material to the
laboratory after the specified dates, a
microbiological study of urine may give an
unreliable result.

10. INSTRUCTION FOR USE

I. Preparing finished medium

1. Dissolve 28 g of powder in 1000 ml of distilled
water.
2. Slowly bring to a boil, gently stirring until the
powder is completely dissolved.
3. Sterilize by autoclaving at 121 ° C for 15
minutes (15 psi).
4. Cool to a temperature of 45-50 ° C and
aseptically pour over sterile Petri dishes. If
necessary, additives are added in accordance with
the relevant instructions.

II. Using finished medium (analysis)

1. Work in the laboratory should be conducted
in compliance with SP 1.3.2322-08. At the
beginning of the study, a direct seeding of the
biomaterial on the surface of the Petri dish with
the finished agar, pre-dried in a thermostat to
remove the condensate from the surface of the
agar, should be performed.

Secretion of outer ear

The material from the tampon, obtained when
examined by an ENT doctor, is inoculated on a
Petri dish with agar in a diagonal stroke, then
inoculated with an inoculation loop to produce
isolated colonies.

Urine

The urine is transferred from the container by a
sterile pipette in an amount of 0.1 ml to a petri
dish with agar. The material is carefully ground
over the surface with a Drigalsky spatula.

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BHIWADI- 301019, (Raj.)

Phone No. : 222073, 221337, 225323

CERTIFICATE for in-vitro diagnostic medical device

Nutrient Agar

LOT E3S1F001

Production date 05/2014

No	Parameter	Specification
1	Appearance	
1.1	Consistence	Free flowing powder
1.2	Color	Light yellowish to brownish yellow color
1.3	Absorbability	positive
2	Solubility	28 g/l
3	Transparency and color of the finished media	
3.1	Transparency of the finished media	Clear, with a slight opalescent
3.2	Color of the finished media	Yellow-cream
4	pH of the finished media при +25°C	7.4
5	Toughness jelly of the finished media	810 g / cm2
6	Gelling temperature	34°C
7	Melting temperature jelly of the finished media	82°C
8	Melting duration jelly of the finished media	48 minutes

Growth promotion test Nutrient Agar with ATCC tests

MICROORGANISMS	ATCC	Inoculum density	Growth	Specification
Staphylococcus aureus	25923	10^3 - 10^4	Good-Luxuriant	Yellow or pale-yellow
Escherichia coli	25922	10^3 - 10^4	Good-Luxuriant	Cream

Results after incubation at 35-37°C for 24-48 hours

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LOT E3E1V001

Production date 05/2014

No	Parameter	Specification
1	Appearance	
1.1	Consistence	Free flowing powder
1.2	Color	Light yellowish to brownish yellow color
1.3	Absorbability	positive
2	Solubility	28 g/l
3	Transparency and color of the finished media	
3.1	Transparency of the finished media	Clear, with a slight opalescent
3.2	Color of the finished media	Yellow-cream
4	pH of the finished media при +25°C	7.4
5	Toughness jelly of the finished media	820 g / cm2
6	Gelling temperature	34°C
7	Melting temperature jelly of the finished media	83°C
8	Melting duration jelly of the finished media	48 minutes

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LOT V3E7B001

Production date 05/2014

No	Parameter	Specification
1	Appearance	
1.1	Consistence	Free flowing powder
1.2	Color	Light yellowish to brownish yellow color
1.3	Absorbability	positive
2	Solubility	28 g/l
3	Transparency and color of the finished media	
3.1	Transparency of the finished media	Clear, with a slight opalescent
3.2	Color of the finished media	Yellow-cream
4	pH of the finished media при +25°C	7.4
5	Toughness jelly of the finished media	820 g / cm ²
6	Gelling temperature	34°C
7	Melting temperature jelly of the finished media	83°C
8	Melting duration jelly of the finished media	48 minutes

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For all questions related to the quality of "Agar
nutrient", it is necessary to contact the official
representative of the manufacturer PROFILAB
LLC (199178, St. Petersburg, Small Ave VO,
house 58, lit. A, Office 19 N Phone 677-52-65,
fax 677-52-66)

Titan Biotech Limited (India)
Official representative in the Russian Federation:
PROFILAB LLC
199178, St. Petersburg, Small Ave VO,
House 58, lit. A, of. 19 N
Tel. 677-52-65, fax 677-52-66
E-mail: info@profilab.spb.ru

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