20

140

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01179820

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Shanghai Lishen Scientific Equipment Co., Ltd. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhou Pingfan

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(Date: Jun. 14, 2022)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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
Revision 2

Director General
Shanghai Lishen Scientific Equipment Co., Ltd


Shen Qinhu

10.06.2022

**HealForce ClassII Biological Safety Cabinets with accessories
produced by Shanghai Lishen Scientific Equipment Co., Ltd, China**

Operations Manual

This manual, including photographs and audio materials, may not be reproduced or distributed wholly or partly without the permission of HEALFORCE.

The hardware operator may copy some sections of this Operations Manual, but these sections shall be intended for internal use, for example for prevention of emergency situations. The provisions of these sections are already clearly marked in this Manual.

HEALFORCE does not incur any liability or obligations for any damages, direct or indirect, caused by not meeting the safety requirements provided in this Manual.

HEALFORCE reserves the right to change this OperationsManual at any time without notice.

Safetyadvice

- * Before using the biological safety cabinet for the first time, please, read this manual carefully!**
- * HFsafe and HFsafeLC(Class II, Type A2)biological safetycabinets shall be used by trained and qualified personnel only.**
- * The hardware maintenance shall be carried out by HEALFORCE or an authorized HEALFORCE agent only.**

Name of the medical device

Biological safety boxes of the second class Heal Force with accessories

Variants of execution:

1. HFSafe - 900 as part of:

- biological safety box;
- stand;
- fluorescent lamp;
- UV lamp;
- instruction manual.

2. HFSafe - 1200 as part of:

- biological safety box;
- stand;
- fluorescent lamp;
- UV lamp;
- instruction manual.

3. HFSafe - 1500 as part of:

- biological safety box;
- stand;
- fluorescent lamp;
- UV lamp;
- instruction manual.

4. HFSafe - 900 LC as part of:

- biological safety box;
- stand;
- fluorescent lamp;
- UV lamp;
- instruction manual.

5. HFSafe - 1200LC as part of:

- biological safety box;
- stand;
- fluorescent lamp;
- UV lamp;
- operating manual.

II. Accessories:

1. Pedal switch.

2. Wheels for the stand – 4 pcs.

3. ULPA filter.

4. HEPA filter.

5. Valves - 3 pcs.

6. Sensors - 2 pcs.

7. Pipeline system.

Introduction

First of all, we thank you sincerely for purchasing our product.

To enable you to use the biological safety cabinet correctly, please, read the operations manual carefully. Otherwise, this may result in equipment damage and even endanger the operator's safety. Please, keep this manual after reading so that it is available for future reference at any time.

Only trained technical personnel authorized, trained and certified by HEALFORCE, are allowed to perform maintenance, adjustment and repair of the biological safety cabinets.

★ Devicedescription

Biological safety cabinets provide both reliable protection for the operator and environment, as well as guaranteed protection of the device against all contamination. Biological safety cabinets are Class II laminar downflow air cabinets used to protect personal, device and environmental safety.

The glass front door of the cabinet can be manually slid freely without direct motor drive, sliding is induced by a counterweight. The high-efficiency particulate air filter (HEPA) achieves unsurpassed air purification rates of 99.999% (for particles of 0.3 μm and larger), and the optional Ultra High Efficiency Filter (ULPA) achieves unsurpassed air purification rates of 99.9995% (for particles of 0.12 μm and more), as a result, the cleanliness of the working space exceeds Class 100.

Class II biological safety cabinets are designed to ensure a clean working environment, while the exhaust air is released through a special filter, automatically ensuring protection to people and the environment. Biological safety cabinets can be used in biological laboratories, medical laboratories, pharmaceutical industry and other related branches, designed to improve process conditions and ensure the operator's health advocacy.

The horizontal inlet air flow that occurs in the front working part is intended for protection; the vertical laminar air flow, after highly efficient filtration through HEPA filters, protects the device against possible contaminations. The exhaust air is filtered through the exhaust air filter, which degree of purification is equivalent to the main filter to protect the environment. If you are not particularly familiar with the operation of the biological safety cabinet, go to Section Operation and Safety Instructions.

The HFsafe and HFsafeLC Type A2 Class II downflow cabinets comply with the US National Standard for Biological Safety Cabinets NSF/ANSI49:2005, European Standard for Biological Safety Cabinet EN12469:2000, and are designed to handle biohazard level 1, 2 and 3 agents.

ATTENTION! If virulent chemical substance or radioactive materials are part of the microbiological study, cabinets specially designed for handling such specimens should be used.

HFsafe and HFsafeLC biological safety cabinets are Class I devices. The devices comply with the following standards:

1. ISO9001/ISO13485
2. GB4793.1
3. GB/T14710
4. YY0569-2005
5. YZB/TEUC021-2005

Interpretation of the model names:

HFsafeandHFsafeLC...biological safety cabinet

900working zonewidththis**900** mm

1200working zonewidththis**1200**mm

1500working zonewidththis**1500**mm

★About this Manual

Safety in many respects depends on correct operation, this manual helps you understand the right way to install, use and maintain your biological safety cabinet. The Manual also includes a description of safety precautions and other important information in order to avoid the hardware damage or your health hazard, please, read this manual before use.

Introduction: Briefly describes the biological safety cabinet, the content and structure of the manual, meaning of the designations in the manual.

Application: Describes the areas in which the biological safety cabinet is to be used and the actions to be taken before use.

Operation n safety precautions: Briefly explains the operation of biological safety cabinets and safetyprinciples.

Device description: Diagrams of the HFsafe and HFsafe LC biological safety cabinet and its components are given.

Installation:Detailed explanation of how to properly transport and install the device, including how to choose the right installation location, and how to install the cabinet components and additional accessories.

Operation: This section helps you to properly use the HFsafe and HFsafe LC biological safety cabinet for efficient safety assurance.

Care and Maintenance: This section recommends a maintenance plan, detailing safetyprecautions during operation and maintenance. Please read this section carefully before carrying out repair and maintenance.

Troubleshooting: This section contains a list of malfunctions that may occur during operation, possible causes of malfunctions and their repair.

Appendix A: Contains a list of components for the HFsafe and HFsafeLC biological safety cabinets.

Appendix B: Contains a list of all additional accessories for the HFsafe and HFsafeLC biological safety cabinets.

Appendix C: Contains electrical ratings and other technical data for the HFsafe and HFsafeLC cabinets.

The last part of the Manual is a warranty statement.

★The symbols used in the Manual are described below:

	Note Symbol: This symbol warns the user that this part is important. Information that may help you gain better insight into the device.
	Biohazard symbol. This symbol warns the user about the presence of a biological hazard. Please be sure to read and understand, otherwise it may lead to serious consequences: illness and even death.
	Operation limitation symbol: All activities marked with this symbol may be performed by professionally trained personnel only. Failure to do so may endanger the user's safety or damage the biological safety cabinet.
	Warning symbol: This symbol highlights important information, please read carefully and observe the warning. Failure to comply may result in serious damage. The warranty does not cover damages due to this cause. HealForce bears no responsibility for any damages caused.

Application

★Scope

Biological safety cabinets are widely used in microbiology, biomedicine, genetic research, biotechnology and other areas. It is the main safety device in the laboratory. Its main purpose is to protect the product, operator and environment against biological hazards that may arise during microbiological work.

The HFsafe and HFsafeLC biological safety cabinet shall be connected to exhaust ventilation when working with samples containing volatile toxic substances or radioactive materials.

★Preparation for installation

The biological safety cabinet provides personal, product and environmental protection during microbiological work. It is necessary to meet the following requirements for:

- Installation location
- Power supply network at the installation location
- Ventilation at the installation location
- Service lines
- Space for maintenance of the biological safety cabinet.

The laboratory area of practice determines requirements for the installation location of the biological safety cabinet. When placing the biological safety cabinet, leave space for service maintenance, convenient routing of the electrical network and service lines. The cabinet shall be placed as far as possible from places of frequent transfer of the personnel and air movement (for example, doorways, places near airshafts, blowers, etc.). More information can be found in Section "Installation".

Operation and safety precautions.

All biological safety cabinets operate according to the following principles:

- Air is filtered by high/extra high particulate air filters (HEPA/ULPA).
- By controlling air flows, a negative pressure area is formed around the working zone and a falling laminar air flow is formed on the working surface.


At the same time, for reasons of safety, the biological safety cabinet shall operate under conditions as recommended by the manufacturer. The operator shall also meet strict operating requirements. HealForce does not bear any responsibility for damage caused by non-compliance with safety precautions.

★HEPA Filter (ULPA - optional)

The high efficiency particulate air filter (HEPA) eliminates more than 99,999% of particles $\geq 0.3\mu\text{m}$ in diameter. Optional Ultra Low Particulate Air (ULPA) filter provides 99.995% purification efficiency for $0.12\mu\text{m}$ particles.

Figure 1

HEPA filter circuit

The filters used in these cabinets are based on a fiberglass cloth. This type of fabric is obtained by binding fibers with a hydrophobic adhesive agent. The filter material is accordion-folded to increase the surface area. In order to prevent the filter damage and folding, and ensure a proper airflow, the filter is installed in an aluminum alloy frame. Rough handling may damage the filter attachment points or create holes in the filter fabric, therefore, for safety reasons, the filter should be exposed to leak checks after installation or relocation.

HEPA filters and airflow ensure operator, product and environmental safety. Therefore, the filter integrity and proper operation are very important. The filter material is very fragile. Never touch it. If there is a suspicion that the filter is damaged, do not use the cabinet. Only professionals shall check filters for integrity before using the cabinet.


HEPA/ULPA filters are only **suitable** for particle filtration. Gas can successfully pass through the filter.

★Airflowsarrangement and control

The air flow pattern in HFsafe and HFsafeLC Class II Type A2 downflow biological safety cabinet consists of a blower, downflow, horizontal inlet airflow, exhaust airflow. In a biological safety cabinetan air flow induced by the blower creates an area where positive pressure is surrounded by negative pressure. The air mixes with the polluted air and 70% of the air flows down through HEPA or ULPA filter, resulting in a downward laminar flow of clean circulating air. Thirty percent of the air passes through exhaust HEPA or ULPA filter to the outside.

Meanwhile, the inside air enters the cabinet through the opening on the front part of the cabinet, forming a horizontal inlet airflow. An air curtain formed from the downflow prevents the inside air from reaching the working surface. The air is directed to the negative pressure area and further to the blower through special openings at the front of the working surface. The air flow diagram is shown in Figure 2.

Figure 2.Diagram of air flowsin the biologicalsafetycabinets.


Laminar downflow

The working zone of HFsafe and HFsafeLC Class II Type A2 downflow biological safety cabinet forms a laminar downflow. In a laminar flow particles move in the same direction at the same velocity. After being filtered through HEPA or ULPA filter, clean air enters the working zone forming a laminar airflow (see Figure 3). The flow captures all suspended particles formed on the working surface and mixes with the horizontal flow coming from the room. One part of the mixed air is filtered through the main HEPA or ULPA filter and again forms a laminar flow, the other part through the exhaust HEPA or ULPA filter is directed out of the device.

In order to determine if a correct laminar downflow is formed, it is necessary to perform tests at several working surface points. The deviation of the average downflow velocity at each point shall be less than $\pm 20\%$ of the average flow velocity.

Figure 3

Laminar downflow pattern for the biological safety cabinet


Horizontal inlet airflow and exhaust airflow

Room air through the front operating glass enters the HFsafe and HFsafe LC Class II Type A2 downflow biological safety cabinet. The inlet air forms a horizontal flow. The mixed air is exhausted through HEPA filter at the top of the cabinet (see Figure 4). The volumes of inlet air and exhaust air are always equal. The horizontal inlet flow comes into the negative pressure area through openings in the front part of the working zone protecting the product in the working zone. At the same time, it prevents particles from flying out of the workingzone, effectively protecting the personnel. Before release HEPA or ULPA filtering precludes the environmental pollution. The average inlet air velocity shall reach a certain value for efficient protectionof the personnel.

Figure 4

Inlet air and exhaustflow patternforbiological safety cabinet

Blower

The blower plays a key role in arranging and controlling HFsafe and HFsafe LC Class II Type A2 biological safety cabinet airflows. The blower draws air into the cabinet through the front operating glass and creates a negative pressure. The blower sends air in two directions. One part is sent through the HEPA filter to the working zone, where a laminar downflow is formed, the other part is released through the HEPA filter from the cabinet. To assure safety throughout the cabinet life, the volume of pumped air shall not vary by more than 10% with a pressure fluctuation of no more than 50%.

★Safety precautions

※Safety check

1. The average downflow velocity in the working zone reaches 0.35 m/s;
2. The filters shall meet integrity requirements;
3. The circulation and exhaust system shall not leak.

1. The safety check shall be performed when the cabinet is turned on for the first time, when the cabinet is turned on again after a long stop in operation and after moving the device!
2. In case of a malfunction that indicates the need for a safety check, the safety shall be checked after the malfunction is corrected!
3. After the malfunction is corrected, and leaks are suspected.
4. Safety checks shall be performed after replacing the filter!
5. All local safety requirements shall be met. The filter integrity and flow velocity shall be tested by the local sanitary epidemiological station or the Center for Disease Control and Prevention or other authorized body to ensure safety of the biological safety cabinet!


※ Only persons responsible for operating the equipment are authorized to turn off the biological safety cabinet blower.

※ Only persons who have been trained and instructed are allowed to use this equipment.

※ This device shall be operated from a separate outlet to avoid shutdown if other devices are damaged.

※ Do not operate the equipment when an alarm is triggered, this may result in material loss, contamination and danger to the operator.

※ Used or contaminated filters shall be placed in a sealed bag, have clear labels in a prominent place stating that the contents are dangerous and may contain dangerous microorganisms. Dispose the filters in a proper way.

The height of the open glass door shall not exceed a safe height, otherwise it will endanger the operator. Do not block the emergency air intake port.

Make sure that the rated voltage and frequency of the mains are consistent with the cabinet input power parameters in order to avoid damage to electrical components and inflammations.

Some components inside the biological safety cabinet may become contaminated during operation. Only specialists can decontaminate the cabinets.

Avoid using open flame in the biological safety cabinet. Open flame can disturb the laminar flow. Long-lasting heating may damage the HEPA or ULPA filter. High pressure gas cannot be used in the biological safety cabinet. Therefore, prior to operation the use of equipment shall be evaluated and tested.

It is recommended to use special infrared light sources.

If objects are to be moved or taken out during operation, try to slow down the operating arm to avoid an impact on air motion, inability to do that may result in emission of polluted air. Never place any object on the grids or block the openings, inability to do that will disturb the airflow and cause biological contamination.

Connect the device only to a grounded outlet. Make sure that the device grounding terminal is connected to the outlet ground. The power supply shall be equipped with current leakage protection. Power down when replacing the fuses. Make sure that the new fuse matches the previously installed one.

Filters installed in the biological safety cabinet will gradually accumulate particles from the room air and working surface. The accumulation rate depends on the room air purity, operating time of the cabinet and conducted operations. The average life of the HEPA filter before replacement is 5 years.


The filter medium of the HEPA filter is fragile - do not touch! Avoid situations in which the filter may be damaged during operation or installation. If you suspect that the filter has been damaged, do not proceed using the cabinet.

The cabinet shall be disinfected in the following cases:

- A. Before performing repair operations in the contaminated area.
- B. Before replacing the HEPA/ULPA filter.
- C. Before the functional test in the contaminated area.
- D. Before reinstalling the biological safety cabinet.

Figure 5. Working zone chart

Areas with dimensions of **900x440mm**, **1200x440mm**, **1500x440mm** on the worktop with no holes are safe working zones, the area outside the safe working zones is not protected. When using the device, you shall take care, not place soft, thin objects (such as soft cloths, etc.) on the table, as they may be sucked into the negative pressure channel and blower, or block the grid holes, which will affect the normal air motion.

When the life of the high-efficiency filter expires, the device emits audible and light alarms, the operator shall contact our service department to replace the filter. Long-term use of the filter will affect the safety records of the device. Periodic check of the filter specifications is very important!


Device description

★Design

H Safe Class II Type A2 biological safety cabinets consist of two major components: the cabinet itself and support stand (See Figure 6)

1. The cabinet housing is made of epoxy-polyester coated steel sheet. The working zone chamber is made of high quality stainless steel.
2. The support stand is made of epoxy-polyester coated steel tubes.
3. Negative pressure air ducts are located on the back and on each side. Each air duct effectively reduces pressure by isolating the working zone from the environment, creating a negative pressure area around the working zone. This ensures no leakage even if the seals are damaged.
4. The housing is equipped with a blower, main and exhaust filter, two waterproof outlets and an UV lamp (optional).

Figure 6. Housing and support stand of the HSafe biological safe cabinets.


5. The holes depicted on the rear wall are used to install liquid and gas delivery cocks(option).
6. There are several bolts on the top cover, they serve to secure the adapter (optional) to the external ventilation.
7. The worktop holes are used to intake the inlet air and polluted air from the operating surface.

Use the drain cock of the tray when cleaning only.

Figure 7. Biological safety cabinet block diagram

1. Sealed blower compartment panel 2. Surface-mounted front panel 3. Fluorescent lamp 4. UV lamp 5. Waterproof electrical outlets 6. Exhaust filter 7. Housing 8. Blower 9. Counterweight. 10. Main filter and deflector 11. Thermoanemometer 12. Working surface 13. Support stand 14. Control console 15. Side walls 16. Faceplate 17. Service ports.

★ Air filtration system

The air filtration system is the main system ensuring a correct functionality of the biological safety cabinet. It consists of a blower, air ducts, main and exhaust filters. The air flow arrangement pattern is depicted in Figure 2.

The main function of the air filtration system is to continuously supply clean air to the working zone at a constant speed of $0.35 \pm 0.025 \text{ m/s}$ with a deviation of $\leq 20\%$. The air purity in the working zone reaches Class 100. In addition, the exhaust airflow is purified, thereby protecting the environment.

★ Front glass

The front glass is raised to an operating height and moved down along the rails using a counterweight. To ensure that the correct operating height is set, the cabinet is equipped with a limiter with an indicator and alarm signal. In addition, the front glass can be equipped with fall protection.

The maximum lifting height of the front glass is 400 mm, while it is possible to install necessary tools in the working zone. But for the time, the cabinet is not protected, the system will emit an alarm signal until a safe height is restored.

Protect the glass against shocks. The glass damage causes a danger to the personnel.

★ Lighting

To create comfortable working conditions, the cabinets are equipped with two fluorescent lamps. They are installed in the front part of the working zone on the front panel. The following lamps are installed: 30W (HFsafe and HFsafeLC-900), 36W (HFsafe and HFsafeLC-1200), 58W (HFsafe-1500)).

★ Bactericidal UV lamp (optional)

HFsafe and HFsafe LC Class II A2 downflow biosafety cabinet accommodates 18W UV lamps (HFsafe and HFsafeLC-900), 30W UV lamps (HFsafe and HFsafeLC-1200), 36W UV lamps (HFsafe-1500). The lamp is installed in the front top of the working zone so that the UV light illuminates the entire working zone. Software and hardware protection prevents UV radiation from escaping outside the working zone of the biological safety cabinet.

★Control console

Figure 8. Control console diagram

Control console overview.

1. Only persons responsible for operating the biological safety cabinet can power down using the on/off button.
2. Button **MODE** is used to navigate between functions in the menu.
3. Button **SET** is used to save the parameters values.
4. Button **LAMP** is used to switch on/off the fluorescent lamps. If the lamps are on, the indicator on the button is illuminated.
5. Button **BLOWER** serves to control the blower. If the blower is on – the indicator is illuminated, when the blower does not operate (for example, with the glass closed) the indicator is flashing. If the blower is turned off, the indicator is illuminated.
6. Button **UV** is used to turn on the UV lamp. Only when the lighting is off and the glass is closed, the UV lamp can turn on. When the lamp is on, the corresponding indicator is illuminated.
7. The **ALARM** indicator is illuminated, when unsafe working conditions are created.
8. The **LCD** screen displays the setup menu or alarm signal.
9. The filter status indicator informs about the filter health and system safety.

Figure 8. Screen menu.

1. **Downflow:** Displays airflow rate in m/s. The downflow rate can be displayed in m³/h (changed by user). The value will flash if the current flow velocity does not comply with the set value with a deviation of $\pm 10\%$. The relevant information will be displayed on the System Status screen.
2. **Exhaust flow:** Displays horizontal inlet flow velocity (complies with the exhaust flow) in m/s or m³/s. The displayed value will flash, if it does not comply with the set value with a deviation of $\pm 20\%$. The relevant information will be displayed on the System Status screen.
3. **Temperature:** Displays temperature in the working zone.
4. **Total time the filter has been used:** Displays an overall time the filter has been used.
5. **Circulation mode:** is used to display in turn downflow velocity, exhaust flow velocity, temperature, and total time the filter has been used.
6. **Menu:** is used to enter the setup menu.
7. **Filter status indicator** displays the filter readiness. Flashes when the filter needs to be replaced. The relevant information will be displayed in the System Status field.

Control console handling.

※ When the cabinet is turned on, the screen first displays the software version, See Figure 10. After a while, the main menu is loaded.

Figure 10. Software version.

VEr1.8

The main screen displays the vertical flow velocity, exhaust flow velocity, temperature and total time the filter has been used.

1. **Vertical airflow.**
The vertical airflow velocity is displayed in m/s (Figure 11) or can be presented in the form of airflow rate in m³/h.

Figure 11. Displaying the vertical airflow velocity.

2. **Exhaust airflow (inlet air velocity)**
The exhaust airflow velocity is displayed in m/s (Figure 12) and also can be presented in the form of air rate in m³/h.

Figure 12. Displaying the exhaust airflow velocity.

3. The temperature in the working zone is displayed in °C (Figure 13).

Figure 13. Temperature display.


4. Total time the filter has been use.
The total time during which the filter has been used is calculated and displayed in hours (Figure 14.)

Figure 14. Displaying the total time the filter has been used.

5. Setup menu.
The BLOWER or LAMP buttons are used to scroll through the list of parameters. Use the SET button to enter the editing mode of the selected parameter. The MODE button returns to the main menu.


Figure 15. Setup men.

5.1. Parameter 01: Password

Default password: 00 (if the password is 00, the system has no protection). (Figure 16)

Figure 16. Password menu.

The user can increase or decrease the password value with BLOWER and LAMP keys. Save the changed value by pressing the SET button and press MODE to return to the main menu.

If the password has been changed, the password shall be entered before each menu entry, or 2 minutes after the last action.

5.2. Parameter 02: Airflow velocity display function.

By default, the airflow velocity is displayed in m/s (Figure 17). With the BLOWER or LAMP buttons, the velocity can be changed to airflow rate in m³/h (Figure 18). Use the SET button to save the new parameter value, and the MODE button to return to the main menu.

Figure 17. Displaying the airflow velocity.

Figure 18. Displaying the airflow rate.

★ Airflow velocity sensor system.

The system consists of a downflow air velocity sensor, an exhaust air velocity sensor, a control system temperature sensor. The temperature sensor measures the downflow temperature. At this time, the downflow velocity sensor measures the velocity, which is displayed on the screen after being calibrated by the measured temperature. At the same time, the data is sent to

the control system, which adjusts the blowers so that the downflow velocity is constant at 0.35 m/s. If the downflow velocity deviates by more than 10%, an alarm is issued. The exhaust air velocity sensor measures the exhaust airflow velocity, if the velocity goes beyond the safety limit, an alarm is sent.

★ **Waterproof electrical outlets.**

There are 2 waterproof electrical outlets in the working zone of the biological safety cabinet.

- Voltage is 220V 50Hz. Maximum current is 5A!
- Water protection is ensured only when the lid is lowered!

★ **Fuses.**

The following fuses are installed in the HSafe Class II Biological Safety Cabinet:

F10A Ø5x20mm – 2 pcs. are mounted in the middle inside the front panel - these are fuses for connection to the mains

F5A Ø5x20mm – 2 pcs. are installed on the right inside the front panel – these are fuses of electrical outlets.

★ **Monitoring and alarm system.**

The system monitors downflow and exhaust flow, front glass position, HEPA filter life and other parameters. If A parameter deviates from the safe values, the system sends an audible warning signal, and the screen displays information on the activated alarm.

The downflow velocity is continuously measured and displayed in real time on the screen. If the flow velocity deviates by more than 10% from the set value, the system automatically sends an audible warning signal, the red indicator lights up, the velocity value on the screen starts flashing. When coming back to safe values, the alarm is turned off.

Do not continue work if the downflow velocity is not correct. Take measures to prevent infection, find out the cause.

The exhaust flow velocity is also measured continuously and displayed on the screen. If the flow velocity deviates by more than 10% as compared to the set value, the system automatically sends an audible warning signal, the red indicator lights up, the velocity value on the screen starts flashing. When coming back to safe values, the alarm is turned off.

- Ensure unobstructed air outlet, do not cover the emergency air intake port.

The position of the front glass is monitored by several limit sensors. If the glass is not in a safe position and is not closed, an alarm signal is issued. Observe the correct height of the front glass during the cabinet operation.

During the operation of the biological safety cabinet, the HEPA filter will slowly accumulate dust and microorganisms. Therefore, the filter resistance will increase. When it gets high enough, it should be replaced with a new one. The control console screen displays the total time of using the filter and the remaining filter life. When the service life expires, an alarm signal is issued. It is recommended not to use the cabinet until the filter is replaced.

To use the filters for a longer time, the cabinet should be installed in rooms with air quality higher than Class 10000. Operation in polluted air will result in frequent filter replacements.

When replacing HEPA filters, the biological safety cabinet shall be open. Therefore, before replacing the filters, it is necessary to carry out a complete disinfection of the cabinet.

Installation.

★ Shipping and transportation

HFsafe and HFsafe LC cabinets can be transported at a temperature of $-20 \sim +50^{\circ}\text{C}$, humidity $\leq 90\%$. The HFsafe and HFsafe LC Class II Type A2 downflow biological safety cabinet consists of a housing and stand. Both parts shall be packed and transported separately. Do not tilt or turn over the cabinet during transportation. Transport the cabinet with care to avoid injuries. Never place the cabinet on the edge, hit, or expose to precipitations or direct sunlight.

★ Installation location


1. The operating temperature is $+10 \sim 30^{\circ}\text{C}$, relative humidity is $< 75\%$.
2. The device shall be installed on a secure, even, horizontal surface.
3. The cabinet should be installed near an electrical outlet. If there is no connection to the external ventilation system, the distance to the ceiling shall be at least 200 mm, the distance of both sides to the wall shall be at least 200mm, the distance of the rear part of the cabinet to the wall shall be at least 300mm.
4. When installed at high altitudes, the airflow rate shall be calibrated.

1. The laboratory in which the cabinet is installed shall be clean, at least isolated from external contamination. It is recommended to install the cabinet in laboratories of a certain cleanliness to ensure maximum performance of the equipment and extend the filter life!

2. The cabinet should be installed away from areas of frequent transfers of the personnel, in front of or near the glass and openings! Figure 19 shows an example of the biological safety cabinet locations, where A, D, E are relatively good locations.
3. The laboratory shall be equipped with ventilation equipment, if the tests carried out can lead to contamination, then the exhaust ventilation device shall be equipped with a proper air filter!

Figure 19. Examples of biological safety cabinet locations.


★ Installation process

1. Carefully unpack the equipment. Check the completeness according to the packing list.
2. Assemble the support stand as shown in Figure 20.
3. Make sure that the support stand is properly assembled with no gaps or deflection. Then place the cabinet on the stand vertically so that the holes for the fixing bolts in the stand and the bottom of the cabinet line up. Remove the worktop transport locks (save them for future transport). Place a stainless steel worktop in the working zone, with the side with a good many vent holes facing the front glass.
4. Remove the counterweight transport lock, place it in the installation location (See Figure 20)
5. Connect the power cable to a fixed electrical outlet. Work is completed.

1. Cabinet is 220kg(HFsafe and HFsafeLC-900), 260kg(HFsafe and HFsafeLC-1200), 300kg (HFsafe-1500). When installing, please pay attention to safety!
2. Be careful when installing, do not damage the seals, damage to the seal may lead to the leakage risk!
3. The cabinet should be installed close to an external fixed electrical outlet, and the fuse of the electrical outlet should not be less than 16A!
4. Check the vent valve seal in the bottom of the cabinet, close the valve!

Figure 20. Biological safety cabinet block diagram

1. Left part of the support stand. 2. Rear part of the support stand. 3. Right part of the support stand. 4. Transportation position of the counterweight fastener. 5. Counterweight fastener location.

★ Installation of gas cocks and other service cocks.

After the cabinet is correctly installed, remove the plugs of three service holes on the rear wall of the cabinet, now service cocks can be installed.


According to Figure 21, place the cock seal washer into the service port, place the cock fitting into the washer, screw the cock to the fitting.

It is possible to install special cocks and service cocks in other places of the working zone at the customer's request.

When connecting liquid, gas or vacuum cocks, pay special attention to tightness of the connections and lines to exclude possible leakages.

Figure 21. Block diagram of additional cabinet connections

A.: Installation of service cocks (vacuum, gas, liquid)

1. rubber bush
2. nut
3. hose pipe
4. clamp
5. hose pipe body
6. external pressure ring
7. rubber sealing washer
8. hose pipe nut
9. latex pipe.

B: Installation of liquefied gas cock.

1. cock body
2. latex pipe
3. valve connector
4. cock washer
5. clamp

C: Connection to the external ventilation

1. clamp
2. adapter to the external ventilation.

★ Installation of adapters in the external ventilation

1. Place the adapter on the emergency air intake port and attach with screws (Figure 21).
2. After installing the adapter, connect the connecting pipe (prepared by the user) and secure with a clamp. Check the interface integrity.
3. The recommended connecting pipe diameter is 200mm. When using pipes of a different diameter, cabinet air flows shall be calibrated,

1. After connecting to the external ventilation, air flows in the biological safety cabinet shall be calibrated!
2. It is recommended to install a check valve on the air duct!
3. When connected to the ventilation system, ensure a sufficient space for installation of air ducts on the biological safety cabinet. When installing, avoid crushing of the air ducts.

★ Startup and usage of the biological safety cabinet

1. Insert the plug into a fixed electrical outlet.
2. Put on a clean laboratory coat and wash your hands. Thoroughly wipe the working surface with 70% alcohol or other disinfectant.
3. Press the "POWER" button (If password is set, enter the correct password to access the system).
4. Check the parameter values on the control console screen.

5. Close the front glass and turn on the UV sterilization lamp. Turn off the UV lamp after sterilization (if sterilization is not required, do not turn on the UV lamp).
6. Switch the safety cabinet in an operating mode. Required parameter values
DOOR: SAFE(the front glass at a safe height); BLOWER:ON (the blower is on); LAMP: ON (the lamp is on); UV: OFF (the UV lamp is off)
7. After the cabinet has operated normally – SystemSafe status, a safety test is performed.
8. If the cabinet has passed the safety test, it is ready for use. Make an appropriate entry in the operation log.
9. The working zone of the safety cabinet shall not be overloaded, all items and tools shall be put in their places before starting operations in the cabinet. Avoid items movements out of or into the working zone during operation - this affects airflows. **All items shall be placed on the working surface in the area with no holes.**
10. Upon completion of the work, all materials and other items used in the process of work shall be disinfected before being removed from the biological safety cabinet. Wipe all workingzone surfaces with 70% alcohol.
11. Do not turn off the cabinet for some time upon completion of works (depending on the material hazard). Close the front glass and turn on the UV sterilization lamp.
12. Turn off the blowers to stop air circulation.
13. Set a password to activate the security (if required).
14. Turn off the fluorescent lamp, close the front glass.
15. Press the **POWER** button to turn off the cabinet.
16. Record all actions in the "operation log"

- 1.Point 7 applies only to cases where safety tests shall be performed!
- 2.If you set a password, remember it correctly.

3. If the alarm indicating the HEPA filter efficiency loss is activated, immediately stop the operation and check all other cabinet parameters. Inform the local sanitary and epidemiological stations, the Centers for Disease Control and Prevention, or other authorized organization for safety verification. If the filter fails during the test, contact the service department to replace the air filter.

If there is a need to bring something in or out of the working zone, move your hands smoothly and as slowly as possible. This will enable to avoid the disturbance of airflows in the cabinet.

★Other control functions

I. Sliding the front glass door

The sliding front glass door is freely lifted and lowered down by hand. If the glass door position is above the safe level during the lift, an alarm will be activated, and the screen will display the front glass overheight warning as shown in Figure 22.

Figure 22. Message:exceeding the rise height of the front glass.

If the glass door position is below the safe position, an alarm signal will be activated, and simultaneously a low level glass message will appear on the screen, as shown in Figure 23.

Figure 23. Message: low level of the front glass.

If the front glass is in a safe position, the display window displays the main screen.

If the front glass is closed, a corresponding message will appear on the screen (Figure 24). At this, the blower will automatically turn off. The UV lamp can be turned on after the fluorescent lamp is turned off.

Figure 24. Message: front glass closed

1. Make sure the front glass door is in a safe position, otherwise, it will unpredictably affect the safety of the equipment.
2. Lower down the front glass door completely before turning off upon completion of the work to prevent contamination of the working zone.

If you close the front glass door while the blower is running, the blower will turn off automatically and the blower indicator will start flashing. When the front glass opens it will automatically start the blower, and the blower indicator will be on continuously.

II. Blower control button **BLOWER**

Controls the blower start and stop. If the blower is off (blower indicator is off), press the **BLOWER** button to start the blower. After starting the blower, if the front glass door is closed, the blower will not start, while the blower indicator flashes, opening the front glass door will turn on the blower, the blower indicator will be on continuously. If the blower is on (blower indicator is on or flashing), press and hold the **BLOWER** button for 5 seconds to turn off the blower. The screen displays the blower off, see Figure 25.

Figure 25 Message: the blower turned off

1. When the blower is turned on, the system will perform mandatory self-cleaning for at least 2 minutes, the screen will display "CLEAN"
2. When the blower is turned off, the blower indicator will go out, but the blower will also perform a mandatory self-cleaning for 2 minutes.
3. Do not turn off the blower during operation. This is a critical factor ensuring the efficiency of the biological safety cabinet.

III. Lamp control button **LAMP**

The lamp is designed to turn on and off the fluorescent lamp.

When the fluorescent lamp is on, the fluorescent lamp indicator is constantly on. When the lamp is turned off, the indicator goes out.


The lamp ensures comfortable illumination of the working zone. Do not turn off the light unless absolutely necessary!

IV. UV lamp control button

It is designed to turn on and off the UV lamp on a regular basis. Each time you press the button, the UV lamp will enter a new operating mode. These modes are repeated in the following sequence: disinfection for 30 minutes, disinfection for 60 minutes, disinfection not limited in time, turning off the UV lamp.

If the fluorescent lamp is on or the front glass door is open, when pressing the UV button you cannot turn on the UV lamp. When the UV lamp lights up, the UV indicator is on.

Figure 26 Setup diagram for UV lamp activation

Display window	Explanation
	disinfection for 30 minutes
	disinfection for 60 minutes
	disinfection not limited in time


Press **MODE** button to exit the UV lamp timer mode.

1. The UV lamp is intended for sterilization, radiation is harmful to human health. By default the UV lamp is off. Do not turn on the UV lamp unless necessary.
2. Before turning on the UV lamp, close the front glass completely. If the operator needs to stay near the cabinet when the UV lamp is on, wear protective clothing and goggles.
3. When disinfecting with the UV lamp, place obvious warning signs outside the location of this equipment!

When the fluorescent lamp is on, or the front glass door is not fully closed, you cannot turn on the UV lamp. In the specified modes, the UV lamp is disabled!

V. **POWER** on/off button is designed to turn the cabinet on/off. If the cabinet is off, press the **POWER** button to turn it on. If the cabinet is running and the user wants to turn off the system, first turn off the blower and close the front glass door, and then press the **POWER** button to turn off the biological safety cabinet. If the glass door is not closed, the display will show the following message (see Figure 27).

Figure 27 Message: the cabinet is turned off.

This method allows users to safely turn off the cabinet. Before turning off, all settings are saved, and the system will use them when it is turned on next time. Use the Power button to turn off the cabinet, and only then you can remove the plug from the outlet.


VII. Blower self-adjustment

The cabinet control system automatically adjusts the blower rpm, depending on the load change to ensure a constant laminar air flow velocity. When the front glass is fully closed, the blower turns off to protect equipment and assure safety.

VII. Password protection

This feature prevents an unauthorized access to the cabinet settings and turning on/off, thus increasing system safety.

If the user password is not 00 (by default), then the password protection function is active. Otherwise, the protection is not active. When the password protection function is active, a correct password is required when making any changes to the settings. The screen will display a message, as shown in Figure 28.

Figure 28 Password entry interface

Press the **BLOWER** or **LAMP** buttons and enter the password. Press the **SET** button within 5 seconds to confirm the password.


If the password is incorrect, the user can make settings. Otherwise, the system will prompt you to reenter the password.

Within 2 minutes after entering the correct password, the protection is automatically disabled, allowing to change the settings. After 2 minutes, the protection will be automatically activated, and you will need to enter the password again to make changes to the settings.

To delete the password, enter the function menu, within 2 minutes after entering the correct password, select function 01, the display window will show a message as depicted in Figure 29, press the SET button to confirm.

Figure 29 Password deleting interface

If the user has forgotten the password and cannot operate the biological safety cabinet, it is possible to disable the password. Unplug the power cord, then press and hold the **MODE**, **SET**, **POWER** buttons simultaneously, plug in the power cord again. The audible warning signal will be activated, and the indicator will flash. Hold the buttons in this position for 30 seconds until the alarm indicator turns off and the audible warning signal stops. After that, the password is enabled, and you can freely use the biological safety cabinet.

★ Recommendations for comfortanceuse

Many factors can affect the working efficiency and your health and safety. The following recommendations can help you create a safer, more comfortable, and more efficient operating environment.

Equipment adjustment

The position of the body relative to the working surface and the posture greatly affect the convenience during operations. We recommend you to use an adjustable chair and footrest so that you can choose the most comfortable position.

Don't sit in a fixed position all day long.

Hips, legs, knees

Use an adjustable working chair so that your feet are on the floor or on a comfortable footrest. If you are using a footrest, it should be wide enough to accommodate your feet comfortably in different positions.

Leave enough space under the cabinet to accommodate your legs and knees and allow for comfortable repositioning.

Do not place cabinets or other items under the cabinet that might get in your way.

Hands, wrist

Keep your hands, wrists straight - in a natural position. Do not rest your hands on the working surface while working. Do not rest your palms while working, this can bend your wrists back and put excessive pressure on your joints.

Shoulder, elbow

Adjust the height of your working chair so that your shoulders and elbows can fully relax. Particularly avoid a position in which the shoulder muscles are stretched during a continuous work.

Eyes

Long periods of work easily lead to eye fatigue.

Recommendations:

1. Let your eyes rest at regular intervals.
2. From time to time, look away from the working area, look into the distance to eliminate fatigue.
3. Wipe the front glass regularly, keep it clean.

Recommendations for work

Items that you frequently use shall be placed directly in front of you to avoid frequent head turns.

The items used shall be conveniently located at a minimum distance from you.

When moving items, do it with your whole hand, not just your wrist.

There are many factors which can improve your comfort. Your practice is also important.

When you work for a long time, you need to take a short break, at least every hour. Regular short breaks increase work efficiency as compared to rare long breaks. During the break, stand up and stretch. Stretch your muscles and joints, especially when you often stay in the same position and the same position while working.

Care and maintenance

★ Maintenance plan

The performance capability and safety of the device can only be guaranteed when inspection, maintenance and repair are carried out by qualified personnel.

※ Maintenance intervals

Maintenance shall be performed every 1000 hours of operation or annually and each time at start-up.

Total running time statistics will have a direct impact on maintenance requirements, therefore, we recommend that you record the cabinet operating time in detail for reference and maintenance requests.

※ Recommended maintenance plan

1. Cleaning

Under normal conditions the biological safety cabinet surfaces can be wiped using common detergents. Do not wipe the glass with abrasive products or products that can cause such an effect.

2. Daily and weekly cleaning

- 1) Disinfection and cleaning of the working chamber are carried out.
- 2) The control console is disinfected and cleaned.
- 3) Using a mild detergent or glass cleaner wipe the external surface of the cabinet and the glass surfaces (if you have a UV lamp, it should be also cleaned).
- 4) Check the device parameters according to the manual.
- 5) Write down the work performed in the "operation log".

3. Monthly or more often as needed

- 1) Wipe off the dust accumulating on the external surface.
- 2) Disinfect the inside parts of the equipment.
- 3) Test the equipment for safety during normal operation.
- 4) Write down the work performed in the "operation log".

4. Annually

An integrated safety check of the cabinet is carried out.

- 1) Qualify the cabinet by a certified engineer.
- 2) Check the reliability of the front glass door of the actuator.
- 3) Perform an integrated safety check.
- 4) If a UV lamp is installed, it should be replaced.
- 5) Write down the work performed in the "operation log".

5. Every two years.

Replace the fluorescent lamp

※DOPtesting (“Smoke Test”)

The biological safety cabinet is equipped with DOPport for testing purposes.

The connected pipeline with a smoke generator shall have a positive pressure. Do not turn on the generator until the pipeline is connected. DOP scanning of the HEPA filters shall be performed by trained and authorized HEALFORCE personnel.

DOPtest procedure:

1. Open the cabinet front panel door.
2. Connect the aerosol photometer to the DOP port in front of the HEPA filter and take measurements.
3. After completing the test, close the DOP port and lower down the front panel.

★ Repairandreplacement

1. Maintenance of this equipment is carried out only by specially trained HEALFORCE technical personnel!
2. In the event of a failure that cannot be remedied by the operator, contact the authorized service immediately. Do not repair yourself - it can be hazardous.
3. Do not replace a part without our prior approval. You can order parts from the HEALFORCE technical service department. When doing so, please indicate which model of biological safety cabinet you have purchased and its serial number!

※Replacement of electrical components

- 1) All operations with electrical components shall be carried out by a qualified electrician in a safe environment (in off state).
- 2) Trained maintenance personnel can be responsible for replacing the fluorescent and UV lamp, counterweight and fuses.
- 3) To repair electrical components, the front panel shall be open (See the front panel disassembly)

***Filter replacement**


1. The filter replacement may result in contamination. Be careful when replacing the filter, this is the most critical maintenance procedure for the biological safety cabinet. In addition, you shall be familiar with relevant laws, standards, regulations and design of this biological safety cabinet.
2. Only authorized HEALFORCE trained specialists can do this job, this is the ideal way to minimize the risk of contamination!
3. Use filters certified by HEALFORCE only, other filters are not permitted, as this may lead to unknown risks.

1. Replacement of the high efficiency particulate air filter

- 1) Before replacing the HEPA filter, make sure that the cabinet has been cleaned and disinfected, and all safety-related operations have been completed.
- 2) Check the safety of all filter locations where contamination may occur.
- 3) Seal the hazardous area and working zone against aerosol leaks. Open the filter compartment and spray the fixing detergent over the entire surface of the filter. Carefully remove the filter from the installation location. If necessary, destroy the filter material and pack it in a special sealed container for further disposal. Carefully install a new filter in the reverse order. After sealing the HEPA filter compartment, remove the working zone seal. See Figure 30 for replacement procedure.
- 4) When changing the filters, personnel shall be equipped with protective clothing, gloves and masks.
- 5) The tight pack with the old filter shall be marked as "hazardous waste" and disposed of according to their category.
- 6) After replacement, check the safety of the work area and equipment.

Figure 30. HEPA filter replacement procedure

 - To replace the blower, main and exhaust filters, remove the airtight front panel of the blower compartment.


1. Blower replacement.

1. Before replacing the blower and associated internal components, make sure that the equipment has been cleaned and disinfected and all safety-related work has been completed.
2. Open the sealed front panel.
3. Loosen the clamp on the factory blue air duct of the exhaust filter.
4. Loosen the connecting screws between the blue air duct and the blower.
5. Loosen four M8 hex screws securing the blower (for Model 1800, eight M8 screws).
6. Remove the blower, install a new blower.
7. Seal the cabinet housing.

※Disassembly of the control console (See Figure 31).

- 1) Turn off the equipment.
- 2) Using a crosshead screwdriver remove the fixing screws (left, right sides, on the bottom of the front panel).
- 3) Hold the edges of the control console with both hands, open the control console outward upward until the control console is locked by the locking mechanism, in so doing, you will hear a token sound.
- 4) After that, two persons shall perform the following operation: hold the console with one hand, press the clamping spring of the locking device with the other hand to release the console.
- 5) Slowly lower down the front panel and set the fixing screws back to their original positions.

Figure 31. Front panel disassembly unit.

★ Cleaning and disinfection

To assure safety when handling the biological safety cabinet, it shall be disinfected on a regular basis. Before using the safety cabinet or carrying out any maintenance, disinfect the entire cabinet, mainly by fumigation.

The main fumigation disinfectant is formalin, the neutralizing agent is ammonia. The operation is performed only by a trained specialist.

1. Remove all items from the safety cabinet, place the sterilization equipment on the working surface, and make the appropriate settings. For specific steps, please refer to the Biological Safety Cabinet Sterilization Equipment User Manual developed by our company.
2. Cover, seal the emergency air intake port, close the front glass, seal the cabinet.
3. Proceed with the sterilization operation in accordance with NSF49:2005.

Do not clean and disinfect without informing the service engineer

- ⚠ Disinfectant containing chlorine is corrosive to metal and shall not be used to disinfect the inside or outside parts of the device.
- ⚠ After disinfecting the inside or outside parts of the device with a disinfectant containing chlorine, immediately rinse the disinfected areas with sterile distilled water.

Malfunctions and troubleshooting

The following table lists some of the malfunctions that may occur, possible causes, and remedies. If you cannot solve the problem using the table, please contact our company or our company's local authorized supplier.

Malfunction	Cause	Remedy	Remark
1. The blower does not work, but its indicator lights up	Connection error	Check the blower connection	
	The blower is overheated	Stop the equipment for a while	
	The front glass is closed	Open the front glass	
2. The light or UV lamp is not on	Connection error	Check the connection	
	The counterweight is damaged	Replace the counterweight	
	The lamp is damaged	Replace the lamp	
3. The buzzer sounds continuously, the alarm indicator lights up constantly	Filter fails	Replace the filter	Safety check is required
	EEPROM reading error	Replace the memory module, check the connection	
	The front glass is not in a safe position.	Move to a safe position	
	Sensor fails	Replace the sensor	
	The air duct is inhibited	Clean	

1. If a malfunction occurs, you cannot continue your work. Air circulation shall continue while the cabinet is cleaned. Remove the materials and accessories from the working area and dispose them of properly. Then contact the customer service!
2. Malfunctions and corrective actions carried out shall be recorded in detail in the "operation log" and signed by the person responsible for the device!

If there is a problem, please contact us, we will help you.

Appendix A: Spare parts

HFsafe n HFsafe LC-900 Biological Safety Cabinet

Accessories
Foot - switch
Standwheels
ULPA filter
HEPA filter
UV lamp
Valves
Sensors
Pipeline system

HFsafe n HFsafe LC-1200 Biological Safety Cabinet

Accessories
Foot - switch
Standwheels
ULPA filter
HEPA filter
UV lamp
Valves
Sensors
Pipeline system

HFsafe-1500 Biological Safety Cabinet

Accessories
Foot - switch
Standwheels
ULPA filter
HEPA filter
UV lamp
Valves
Sensors
Pipeline system

Appendix B: Accessories

Accessories
Foot - switch
Standwheels
ULPA filter
HEPA filter
UV lamp
Valves
Sensors
Pipeline system

Appendix C: Performance data

Model	HFsafe and HFsafe LC-900	HFsafe and HFsafe LC- 1200	HFsafe -1500
Dimensions WxHxD, mm	1023x2260x800	1323x2260x800	1623x2260x800
Dimension of the working zone WxHxD, mm	900x678x600	1200x678x600	1500x678x600
Filter	HEPA or ULPA	HEPA or ULPA	HEPA or ULPA
Inlet air flow velocity, m/s	0,53±0,025	0,53±0,025	0,53±0,025
Downflow velocity, m/s	0,35±0,025	0,35±0,025	0,35±0,025
Voltage	220 V 50 Hz	220 V 50 Hz	220 V 50 Hz
Power	800 W	800 W	1200 W
Maximum power	1850 W	1850 W	2200 W
Fuses	2x10A, 2x5A	2x10A, 2x5A	2x10A, 2x5A
Noise level, dB	<65	<65	<65
Weight, NET, kg	220	260	300
Operating temperature	+10~ +30°C	+10~ +30°C	+10~ +30°C
Operating humidity	30-75%	30-75%	30-75%

The above parameters are guaranteed when the operation and maintenance are performed in accordance with this manual, and the original accessories and spare parts are used!

If a user wants to use accessories supplied by our company (for example, footswitch, external ventilation adapters, service cocks, etc.), please clearly indicate the required accessories when ordering the biological safety cabinet. When replacing parts, please contact HealForce or your local supplier and specify the model. Otherwise, we cannot guarantee the correct operation of the cabinet!

Warranty statement

- The warranty period of the device is 12 months from the date of purchase.
- The company bears no warranty obligations for any damage or failure of the device and equipment caused by misuse of equipment.
- After the warranty period, the company provides service support, but at extra charge.
- The company will provide trained and authorized engineers of the company with necessary drawings and technical data.

DISPOSAL AND DESTRUCTION PROCEDURE



The biological safety cabinet shall be disposed off through special companies in accordance with official disposal requirements in your region.

Blower parameters:

Biological safety cabinets	Blower	Diameter(mm)	Power (W)	Intake air volume (m ³ /h)
HFSafe - 900	D4E 225-CC	225	650	2650
HFSafe - 1200				
HFSafe - 900LC				
HFSafe - 1200LC				
HFSafe - 1500	D4E 225-DH	225	1060	3750


Front glass door parameters and specifications:

- the front glass door lifting system makes it possible to fix the glass in two positions, it is raised and lowered down using a system of counterweights
- safe height- up to 200 mm
- counterweight design: the front glass door is equipped with counterweights to facilitate lifting and fixing.

Biological safety cabinets	Overall dimensions (mm) (WxHxD)	Safe height (mm)	Counterweight design
HFSafe - 900	9.2x777x994	200	
HFSafe - 1200	9.2x777x1294	200	
HFSafe - 1500	9.2x777x1594	200	
HFSafe - 900LC	6x670x983	200	
HFSafe - 1200LC	6x670x1283	200	


Diameter of holes for inlet air intake

An opening is provided for the intake of inlet air, with a size of 6x16 mm


Stand dimensions

Biological safety cabinets	Stand dimensions			
	Length (mm)	Width (mm)	Height (mm)	Pipe diameter(mm)
HFSafe - 900	1023	768	610	40x40
HFSafe - 1200	1323	768	610	40x40
HFSafe - 1500	1623	768	610	40x40

Biological safety cabinets	Stand dimensions			
	Length (mm)	Width (mm)	Height (mm)	Pipe diameter(mm)
HFSafe - 900LC	1038	710	660	50x50
HFSafe 1200LC	1338	710	660	50x50

Strength characteristics of additional cabinet connections:



- service cocks a cock made of H62 grade copper (tensile strength is 680MPa, yield strength is 450MPa)
- external ventilation - not applicable

Accessories parameters and specifications:


The air duct design:


- a pressure range in the air ducts is 1000 Pa
- a pressure range in the working zone –atmospheric pressure

Pipelinesystem:

Revision 2

Filters

Biological safety cabinets	Dimensions (mm)	Number of filters (ea.),	Bacteria/microbial permeability	Water-resisting property
HFSafe - 900	550x400x69	1	99.9995%	fiberglass, non-water-resistant
	915x450x69	1		
HFSafe - 1200	720x400x69	1		
	1219x450x69	1		
HFSafe - 1500	890x400x69	1		
	1525x450x69	1		
HFSafe - 900LC	550x400x69	1	99.9995%	
	915x450x69	1		
HFSafe - 1200LC	720x400x69	1		
	1219x450x69	1		

Number and sequence of layers in the filter:

- in the HFSafe- 900 and HFSafe- 900LC cabinets: 223 layers in the main filtering and 145 layers in the secondary filtering;
- in the HFSafe- 1200 and HFSafe- 1200LC cabinets: 298 layers in the main filtering and 175 layers in the secondary filtering;
- in the HFSafe-1500 cabinets: 375 layers in the main filtering and 216 layers in the secondary filtering.

The filters are based on glass fiber (glass paper) in one accordion-fold layer.

Wheel parameters:

- diameter (mm) – 100
- height from floor (mm) – 116
- locking system - the wheels are blocked by pressing the fixing brake.

Parameters and specifications:

- **foot-switch** – water resistance level IPX8 (long-duration immersion in water will have no adverse effect on the device under standard conditions); driving force: 10N~50N

- **Valves** – size G3/8", operating pressure -10 bar,

- **Sensors** – airflow measurement - resolution 0.01 m/sec, range 0~2 m/sec.

Brands of materials used for manufacture of medical devices and accessories:

Steel grade for the cabinet housing

The cabinet housing is made of steel SPCC.

Steel grade for the working zone chamber

The working zone chamber is made of stainless steel SUS304.

Steel grade for the stand

Stand (support) for the cabinet is made of steel Q235.

Epoxy-polyester coating grade

The steel body of the cabinet is coated with RAL9016, manufactured by "Jinhu Color Powder Coating Co., Ltd", China

Grades of materials for filter

The HEPA H14 filter is made of microfine and ultrafine glass fibers (glass paper). The ULPA U17 filter is applicable for a higher degree of purification. It also makes it possible to purify the air from very small particles, including viruses. The filter is made of filtering medium consisting of micro-thin and ultra-thin glass fibers.

Material brands for filters HEPA H14, ULPA U17:

Filter body - aluminum alloy AFF & CAMFIL, manufactured by American Air Filter Company, Inc., the USA, Camfil Gourp, Sweden

The filtering medium is glass fiber (glass paper) of Chuansha brand, manufactured by Chuansha Shanghai Co. Glass Products Co., Ltd, China

Sealant - polyurethane Momenite, manufactured by Momenite Performance Materials Inc, the USA and China.

Grades of materials for the front glass

The material used for the front glass is UV-absorbing tempered glass, laminated, colorless and transparent, at least 5 mm thick. The brand of glass used in production is M4, manufactured by Chuansha Shanghai Co. Glass Products Co., Ltd., China.

Name and brand of materials used

Name of item	Material
Foot – switch	ABS Plastic. Brand: PolylacPA-777C, manufactured by ChiMei, Taiwan/China
Stand wheels	Rubber STR-20, manufactured by SRITRANG GROUP, Thailand
ULPA and HEPA filter housings	Aluminum alloy AFF&CAMFIL, American Air Filter Company, Inc., USA/Camfil Gourp, Sweden. Sealant - polyurethane, Momentive, manufactured by Momentive Performance Materials Inc, the USA and China
ULPA filter	Filter material - Glass fiber (glass paper) Chuansha brand, manufactured by Chuansha Shanghai Co. Glass Products Co., Ltd, China
HEPA filter	Filter material - Glass fiber (glass paper) Chuansha, manufactured by Chuansha Shanghai Co. Glass Products Co., Ltd, China
UV lamp	Actinic TL-D manufactured by Philips, China
Valves	Cooper. H62
Sensor	ABS Plastic. Brand: PolylacPA-777C, manufactured by ChiMei, Taiwan, China
Pipeline system	Plastic PVC (polyvinylchloride). Brand: SG-8, manufactured by XINJIANG TIANYE FOREIGN TRADE CO LTD, China
Housing	Sheet steel SPCC. Epoxy-polyester coating RAL9016, manufactured by Jinhu Color Powder Coating Co., Ltd, China
Working zone chamber	Stainless steel SUS304
Front glass	UV-absorbing tempered glass, laminated, colorless and transparent, at least 5 mm thick. Glass M4, manufactured by Chuansha Shanghai Co. Glass Products Co., Ltd., China.
Stand (support) for the cabinet	Steel Q235
Control console	Plastic PVC (polyvinylchloride). Brand: SG-8, manufactured by XINJIANG TIANYE FOREIGN TRADE CO LTD, China
Transport package	Wooden cabinet manufactured by Shanghai Satyr Packing (China)

Purpose of medical device

The biological box is widely used in microbiology, biomedicine, genetic research, biotechnology and other fields. This is the main device for ensuring safety in the laboratory. Its main purpose is to protect the product, the operator and the environment from biological hazards that may arise during microbiological work.

Indications for use, contraindications, possible side effects when using a medical device

Indications for use:

HFSafe Class II biological safety cabinets are indicated for use in microbiology, biomedicine, genetic research, biotechnology and other areas.

It is the main safety device in the laboratory. Its main purpose is to protect the product, operator and environment against biological hazards that may arise during microbiological work.

HFSafe and HFSafeLC biological safety cabinets are applicable to work with microorganisms of pathogenicity (hazard) class I, II, III, samples containing volatile toxic substances or radioactive materials.

Contraindications:

The cabinet is not intended for work with microorganisms of class IV pathogenicity (hazards).

Possible side effects:

If the Operations Manual for medical device HFSafe Class II Biological Safety Cabinet is strictly observed, there are no side effects.

★Quality Policy of the Company:

Steady improvement

Meeting customers' requirements

★After-sales service hotline: 800 -820-6044

Shanghai: Phone: (021) 62527883

Fax: (021) 62130871

★The company has the right to immediately make changes to the equipment specification improving its performance, without prior notice.

★If you have problems with the equipment, please contact the nearest employee of our company, you will be satisfied with the service.

★Shanghai Lishen Scientific Equipment Co., Ltd the owner of HealForce brand

★Manufacturer: Shanghai Lishen Scientific Equipment Co., Ltd

Factory Address: No. 6788, Songze Avenue, Qingpu, 201706, Shanghai

Phone: (021) 59867888

Fax: (021) 59867199


Authorized representative of the medical device manufacturer

Sinopharmtech Limited Liability Company (Sinopharmtech LLC)

Address: 119421, Moscow, Leninsky ave., 105 k. 2, room./com. 20/105

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