UR.-Nr. 1946 / 2017 F

I, the undersigning Notary Public, hereby certify that the signature on the first page of the attached "Technical File" and the signature on the first page of the attached "Instruction for Use" is the true signature, which was signed before me, of

Dr. Ingo Heschel, born on 6th February, 1961, business address at 52134 Herzogenrath/Germany, Kaiserstraße 100,

acting not in his own name, but for the company

Matricel GmbH

with registered office in Herzogenrath/Germany, which is registered with the Companies' Register of the Local Court of Aachen under HRB 8628.

The identity of Dr. Ingo Heschel was proved to me by his passport No. L75TK7CJ2.

Upon today's inspection of the electronic Commercial Register, I certify that

- the company Matricel GmbH is registered with the Local Court of Aachen – Commercial Register – under HRB 8628,
- Dr. Ingo Heschel is entitled to represent the company Matricel GmbH with sole power of representation.

I further certify that the copies of the documents attached in this folder are true copies of documents, which were presented to me today by Dr. Ingo Heschel, who declare that the documents serve for a medical device application in the Russian Federation.

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Left: Inner more dense fibrillar membrane surface; Right: Outer side of the membrane (coarser fibrous structure)

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Raw Material Acquisition

Site of production	Matricel GmbH, Kaiserstrasse 100, 52134 Herzogenrath, Germany
Specification	Compliant with: ISO 22442-1 – Medical devices utilizing animal tissues and their derivatives – Part 1: Application of risk management ISO 22442-2 – Medical devices utilizing animal tissues and their derivatives – Part 2: Controls on sourcing, collection, and handling ISO 22442-3-Medical Devices utilizing animal tissues and their derivatives-Part 3: validation of the elimination and/or inactivation of viruses and transmissible spongiform encephalopathy (TSE) agents
Rationale	Assures safety of raw material
Test Method	Quality Assurance audit of the Device History Record

Material composition

Specification	The DBM device is derived from porcine connective tissues and is composed of a network of highly purified porcine collagen (predominantly Type I and some Type III) and elastin fibers. Calculations show that the percentage of elastin in the harvested membranes averages 15%.
----------------------	---

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process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Do not use if package has been previously opened or damaged.

Packaging

Packaging of the product involves the sealing of preformed blisters with Tyvek lids. The outcome of the sealing process cannot be 100% verified and therefore has been validated. In addition during secondary packaging automated application of the product label on the blisters is performed and has been validated. Qualification of the blister sealer process was performed to establish the optimal sealing conditions and to show that the seal quality is maintained after sterilization under routine EO sterilization (normal condition, 1xEO) and EO re-sterilization (challenge condition, 2xEO). For the blister sealing operation, the predetermined requirements were based on ISO 11607-1 with respect to seal width, seal strength and seal appearance.

Packaging	
Primary, Secondary	Preformed blister - polyethylene terephthalate (PETG), Nelipak, Netherlands Cover Tyvek 1073B by Oliver Tolas, Nelipak, Netherlands
Tertiary	Folding box
Outer	Chipboard box

Dimensions:

Tertiary: 8,2х11х2 (cm)

Appearance of blister: 7,4х10,2х1 (cm)

The inner blister: 5,2х6,8х0,7 (cm)

Weight:

Name	The inner blister	Appearance of blister
15х20 mm	2,09g	6,81 g
25х30 mm	2,23 g	7,00 g
30х40 mm	2,24 g	7,27 g

(Specified max. Values)

All packaging requirements is accepted with ISO 11607

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A yearly evaluation is performed according to the requirements of EN ISO 11135.

Ethylene oxide (EO) sterilization of the DBM device in double blister packaging is performed by the contract sterilizer HA2 Medizintechnik GmbH.

Ethylene Oxide Sterilization of the DBM device

Contract Sterilizer	HA2 Medizintechnik GmbH Am Bahndamm 11 38820 Halberstadt Quality systems certificates: 93/42/EWG, DIN EN ISO 13485:2012+AC:2012, DIN EN ISO 11135:2014
Method of Sterilization	Ethylene oxide
Sterility Assurance Level	10 ⁻⁶
Sterilization Validation Method	ISO 11135-1, Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.
Ethylene oxide residuals (EO and ECH)	ISO 10993-7, Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals

EO = ethylene oxide; ECH = ethylene chlorhydrin

Validation and routine control of sterilization of DBM device by ethylene oxide gas treatment is based on SO 11135-1. This standard describes the requirement to confirm that the product, packaging and labeling materials used are compatible with routine sterilization conditions. In addition, in case it is planned that re-sterilization will be allowed, the influence of re-sterilization conditions on the product, packaging and labeling need to be evaluated.

As part of the sterilization validation activities, experiments were performed to demonstrate compatibility with routine sterilization conditions as well as re-sterilization conditions. Based on the data obtained during this study, it was concluded that the DBM device, the packaging materials and the labeling materials are compatible with routine sterilization conditions as well as re-sterilization conditions.

As described in this overview document, the basic validation of the EO sterilization cycle parameters was performed by HA2 Medizintechnik GmbH (HA2), using a worst case product loading configuration (PVC-loaded dummy pallets). Additional validation activities were performed on the Dental Barrier Membrane, to evaluate different sterile barrier packaging configurations. English translations of the final test reports from HA2 can be provided upon request.

Residual Toxicity

The allowable limits for EO (ethylene oxide) and ECH (ethylene chlorhydrin) are within the published limits of ISO 10993-7:2008. For the purposes of selecting the appropriate allowable limits for residual EO and ECH, the Dental Barrier Membrane is categorized in terms of duration of contact as "permanent contact" (see ISO 10993-7:2008). The allowable limits are set to ensure that:

The average daily dose of EO to patient shall not exceed 0.1 mg/d

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% Other Types Collagen and List of which Types present	Not performed. Matricel Bulk Material Collagen is obtained by removing non-collagenous and non-elastinaceous substances from porcine connective tissue. Therefore the naturally present types of collagen will be present at the naturally occurring ratio in the Bulk Material Collagen. The composition is mainly collagen Type I and low amounts of collagen Type III as for the predicate device Bio-Gide (K042197 and K050446) that uses the same porcine tissue as starting material.
Total DNA (ppm or %)	Not performed. Matricel has selected amino acid analysis to show purity. In addition, the source tissue used (porcine connective) contains only limited cell numbers, making this test method less appropriate.
% native collagen (by trypsin resistance, circular Dichroism)	Not performed. Matricel has selected Differential Scanning Calorimetry as the preferred method to assess collagen denaturation/nativity.
Biodegradable	DBM is intended to be biodegradable. In canine animal studies, an estimated 80-90% of the collagen membrane has degraded within 6 weeks of implantation.

Specifications of the Sterile Finished Product for the Matricel DBM

Inspection Item	Specification	Test Instruction (Matricel Doc. No.)	Test Information
-----------------	---------------	--------------------------------------	------------------

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Remnant production aids:			
pH value	5.0 - 7.4	QS-TI-1016	0.75-1.5 g of the Dental Barrier Membrane is cut into pieces and transferred to a sample tube containing 20 g distilled water and allowed to shake for at least 1 hr (1-4 hrs). The pH of the solution is measured.
Endotoxin level (LAL)	< 20 EU/product	QS-TI-1015	Referenced in – Sterilization and Shelf Life*
Sterility	Sterile based on parametric and bio-indicator release	QS-TI-1045	Referenced in – Sterilization and Shelf Life*
Acetone content	≤ 250 mg/kg	QS-TI-1010	1 g of the Dental Barrier Membrane is cut into pieces and tested for acetone content using Head Space gas chromatography.
Quality and purity:			
Peak denaturation temperature	≥ 49°C	QS-TI-1005	Measure the denaturation temperature specifically of the collagen using Differential Scanning Calorimetry (DSC).
Dimensions/Physical characteristics:			
Color	White/off white	QS-TI-1045	Visual inspection.
Expansion after wetting	≤ 25 %	QS-TI-1053	Expansion is determined by measuring the surface area of the membrane sample before and after wetting with Phosphate Buffered Saline buffer (pH = 7.4 ± 0.1 at 25°C).
Tensile strength	≥ 1.7 N/mm (dry)	QS-TI-1052	Tensile strength is tested using a Z2.5 machine (Zwick/Roell) with load sensor 500 N. Test samples with a width of 5–10 mm

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			and a length of at least 30 mm will be tested. One strip per membrane is measured. The clamping length is 10 mm. The tensile strength determination is performed with a strain rate of 125%/min.
--	--	--	--

Mechanical properties:

A verification study was performed to investigate the mechanical properties of the DBM. The DBM will be manipulated during application and fixed in place according to standard surgical techniques (In dental applications: forceps handling of the dry membrane followed by tissue flap closure with or without additional fixation of the membrane). A minimum requirement for Tensile Strength (TS) of 1.7 N/mm has been defined for the DBM when tested dry. The verification study demonstrated compliance of the device with this requirement.

Tensile strength of membranes as a function of surface weight and sterilization conditions.

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Biocompatibility

The biocompatibility of the DBM device has been tested according to ISO 10993-1, Biological evaluation of medical devices, Part 1: Evaluation and testing was done for an implanted device that contacts tissue/bone for a prolonged period of time. The DBM used in the biocompatibility tests were finished, sterilized devices. All of the biocompatibility testing was performed under GLP (Good Laboratory Practices) conditions at Toxikon Laboratories (United States of

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America) and is referenced in the reports as "collagen membrane – EO sterile". The DBM devices and controls were prepared for testing according to ISO 10993-12:2007 'Biological Evaluation of Medical Devices, Part 12: Sample Preparation and Reference Materials'. Results of the testing support the biocompatibility of the Dental Barrier Membrane for its intended use.

Summary of Biocompatibility Testing

Report No. (Appendix #)	Test	Result
12-2170-G1 (Appendix O1)	Cytotoxicity (ISO 10993-5:2009) (MEM Elution)	Non-cytotoxic
12-2170-G2 (Appendix O2)	Intracutaneous Injection Test (ISO 10993-10:2010)	Non-irritating
12-2170-G3 (Appendix O3)	Kligman Maximization Sensitization (ISO 10993-10:2010)	Non-sensitizing
12-2170-G4 (Appendix O4)	System Injection Test (ISO 10993-11:2006)	Non-toxic
12-2170-G5 (Appendix O5)	Genotoxicity (ISO 10993-3:2003) (Ames test)	Non-mutagenic
12-2170-G6 (Appendix O6)	Genotoxicity (ISO 10993-3:2003) (Chromosome Aberration)	Non-clastogenic
12-2170-G7 (Appendix O7)	Genotoxicity (ISO 10993-3:2003) (Mouse Micronucleus)	Non-clastogenic
12-2170-G8 (Appendix O8)	Intramuscular Implantation: 2 weeks (ISO 10993-6:2007)	Non-toxic
12-2170-G9 (Appendix O9)	Intramuscular Implantation: 4 weeks (ISO 10993-6:2007)	Non-toxic
12-2170-G10 (App O10)	Intramuscular Implantation: 13 weeks (ISO 10993-6:2007)	Non-toxic
12-2170-G11 (Appendix O11)	Intramuscular Implantation: 26 weeks (ISO 10993-6:2007)	Non-toxic
12-2170-G12 (Appendix O12)	Subchronic Systemic Toxicity: 90 Day Systemic Toxicity in Rats via Intramuscular Implantation (ISO 10993-11:2006)	Non-toxic

Cytotoxicity Test (L929 Neutral Red Uptake, 1 Concentration) (Appendix O1):

The biological reactivity of a mammalian cell culture (mouse fibroblast L929) in response to the DBM device (referenced in the report as "Collagen membrane – EO sterile") test article extract was investigated according to ISO 10993-5, 2009. The viability of the L929 cells is determined via metabolic activity as measured by the internalization of Neutral Red in lysosomes of cells. The number of viable cells correlates to the color intensity determined by photometric measurements after extraction of the incorporated Neutral Red.

One lot of the DBM device was used for the study. The test article (60cm²), the positive control (natural rubber, 30cm²) and negative control (high density polyethylene, 30 cm²) were extracted in 10 mL of Minimal Essential Medium with 10% fetal bovine serum (referred to as complete

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MEM) at 37±1°C for 24±2 hrs.

100µL of the test article extract, the control article extracts, and the complete MEM that was subjected to the same time and temperature treatments were used to replace the maintenance medium of L929 cell cultures (grown in 96-well plates) in 6 replicates each (except for the untreated control which is used in 12 replicates). Cultures were incubated at 37±1°C for 24-26 hrs in a humidified atmosphere containing 5±1% CO₂.

After incubation, the culture medium was removed from the plate and the cells were washed with phosphate buffered saline (PBS). Neutral Red medium was added to each plate, and the plates further incubated at 37±1°C for 3±0.2 hrs. The cells were then washed with PBS and an ethanol/acetic acid solution was added to each well to extract the Neutral Red. The optical density of each well was measured at 540nm. A reduction of viability of cells (expressed as viability %) was calculated based on the optical density of the cells exposed to the test article compared to cells exposed to the untreated control complete medium. If viability is reduced to <70% of the blank, the test article is considered to have a cytotoxic potential.

Result: The viability of cells exposed to the neat test article extract was 96%. Based on the criteria of the protocol and ISO 10993-5, 2009, the DBM device is not considered to have cytotoxic potential.

Intracutaneous Injection Test (Appendix O2):

Skin irritation testing after intracutaneous injection of the DBM device test article extract was performed according to ISO 10993-10, 2010. The purpose of the test is to determine if any chemicals that may leach or be extracted from the test article are capable of causing local irritation in the dermal tissues of New Zealand White rabbits (1 male and 2 females were used). The DBM device test article (referenced as “Collagen membrane – EO sterile” in the report) was prepared by extracting 15 cm² of the membrane separately in 20 mL of USP 0.9% saline (NaCl) or cottonseed oil (CSO) at 37±1°C for 72±2 hrs. NaCl and CSO that was subjected to the same time and temperature treatments served as the untreated controls (or blank).

Prior to injection, the animals were examined and the fur clipped. On each side of the animal, 0.2 mL of one of the test extracts was injected into 5 sites along with the control injected into 5 sites giving a total of 20 injection sites per animal. The sites were visibly inspected for erythema and edema immediately following the injection and at 24±2, 48±2 and 72±2 hrs. Observations were scored according to the “Classification System for Scoring Skin Reactions” as described in the attached protocol.

Result: No signs of erythema or edema were noted at any observation periods. The test article sites did not show a significantly greater biological reaction than the control sites. Based on the criteria of the protocol, the DBM device test article when extracted in either USP 0.9% NaCl or CSO is considered a negligible irritant in New Zealand White rabbits.

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Kligman Maximization Sensitization Test (Appendix O3):

Sensitization testing was performed according to ISO 10993-10, 2010. The purpose of the testing was to evaluate the allergenic potential or sensitizing capacity of the DBM device test article extract (referenced in the report as "Collagen membrane – EO sterile") in guinea pigs. The DBM device test article was assessed for sensitization using Hartley guinea pigs (16 males and 19 females). The DBM device was separately extracted in USP 0.9% saline (NaCl) and cottonseed oil (CSO) at 37 ± 1 °C for 72 ± 2 hrs. NaCl and CSO that were subjected to the same time and temperature treatments served as the untreated controls (or blank). 1% DNBC (dissolved in 95% ethanol) served as the positive control.

Ten (10) animals were used to test the test article extracted in NaCl, ten (10) were used to test the test article extracted in CSO, five (5) were used for each negative control and five (5) were used for the DNBC positive control. During each phase, as described in Table, the animals were injected intra-dermally with 0.1 mL of the indicated solution.

Table: Dosing Solutions for Guinea Pig Maximization Sensitization Test

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Phase (Day)	Groups	Treatment
Induction/intradermal Application (Day 0)	<u>Experimental Group</u> (1) 0.1mL FCA + vehicle (2) 0.1mL test article (3) 0.1mL test article 1:1 with FCA <u>Negative Control</u> (1) 0.1mL FCA 1:1 vehicle (2) 0.1mL vehicle (3) 0.1mL vehicle 1:1 with FCA <u>Positive Control</u> (1) 0.1mL FCA 1:1 NaCl (2) 0.1mL DNCB (3) 0.1mL DNCB 1:1 with FCA	3 pairs of intradermal injections were made (6 total injection sites) per animal.
Topical Application (Day 7)	<u>Experimental Group</u> Filter paper saturated with test article extract was wound around the torso of the animal for 48 hrs <u>Negative Control</u> Filter paper saturated with negative control, same protocol as above. <u>Positive Control</u> Filter paper saturated with DNCB, same protocol as above	On Day 6, pretreatment comprising of 10% SDS in petrolatum was given to animals that showed no signs of irritation or corrosion. A 2 X 4cm filter paper was used.
Challenge Application (Day 23)	<u>Experimental Group</u> Filter paper saturated with test article is applied to an unexposed area for 24 hrs <u>Negative Control</u> Filter paper saturated with negative control, same protocol as above <u>Positive Control</u> Filter paper saturated with DNCB, same protocol as above	A 2 X 2cm filter paper was used.

FCA = Freund's Complete Adjuvant
 Experimental = Test article extract in either polar or nonpolar solvent
 Control = Polar or nonpolar solvent alone
 DNCB = 0.1% DNCB in 95% ethanol
 SDS = sodium dodecyl sulfate

Following the challenge application phase, injections sites were observed for erythema and edema at 24, 48 and 72 hrs. Injection sites were scored using the Magnusson and Kligman Scale and the Sensitization Classification.

Result: As defined by the Kligman scoring system, the DBM device test articles extracted in USP 0.9% NaCl or CSO were considered a Grade I reaction and is therefore classified as having weak allergenic potential.

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Systemic Injection Test (Appendix O4):

The systemic injection testing to test for potential irritation effects of the DBM device test article in twenty (20) Albino Swiss mice was performed according to ISO 10993-11, 2006. 15 cm² of the DBM device test article (referenced as "Collagen membrane – EO sterile" in the report) was prepared by extracting it in 20 mL of either USP 0.9% saline (NaCl) or cottonseed oil (CSO) at 37±1°C for 72±2 hrs. NaCl and CSO that was subjected to the same time and temperature treatments served as the untreated controls (or blank). Animals were injected according to Table. Each group consists of five (5) animals.

Table: Dosing Administration for Systemic Injection Test

Group	Dose/kg, route	Injection rate
Experimental Group (test article extracted in NaCl)	50mL, IV	2 mL/minute
Experimental Group (test article extracted in CSO)	50mL, IP	2 mL/minute
Control Group (NaCl)	50mL, IV	2 mL/minute
Control Group (CSO)	50mL, IP	2 mL/minute

IV = intravenous
 IP = intraperitoneal

The animals were observed for clinical and toxicological signs at 4, 24, 48 and 72±2 hrs after injection. During the 72 hr observation period, no animals showed any adverse physical symptoms when compared to the controls.

Results: The DBM device test article extracted in either NaCl or CSO did not induce a significantly greater biological reaction than the control extracts when tested in Swiss Albino mice. Therefore, the results are negative and the test article is considered non-toxic.

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The purpose of this test is to determine the ability of the DBM device test article extract to induce reverse mutations in histidine (his- to his+) and tryptophan (trp- to trp+) genes in four strains of *S. typhimurium* and one strain of *E. coli* respectively. ISO 10993-3, 2003 was followed. The testing was conducted in the presence and absence of an exogenous mammalian metabolic activation system using plate incorporated method of exposure.

The DBM device test article (referenced as “Collagen membrane – EO sterile” in the report) was prepared by extracting 15 cm² of the membrane separately in 5 mL of USP 0.9% saline (NaCl) or cottonseed oil (CSO) at 37±1°C for 72±2 hrs. NaCl and CSO that were subjected to the same time and temperature treatments served as the negative controls.

The bacteria were exposed to the test article extracts with or without a metabolic activation system, allowed to incubate and revertant colonies were counted compared to those of the negative control article cultures. Crystal violet sensitivity and ampicillin resistance tests were used to test requirement of histidine and tryptophan in the bacterial colonies. The test article extracts were applied to each plate at a concentration of 100 µL/plate. Negative and positive controls with or without metabolic activation were prepared according to the protocol. The mean number of revertants per plate was calculated for the test extract, negative and positive controls. Generally, a response that is twice that of the negative control article is considered a mutagenic response.

Results: The DBM device test article extracted in either NaCl or CSO did not induce a significant increase in the number of revertant colonies compared to the corresponding negative controls in both non-activated and activated conditions. Therefore, the test article is considered non-mutagenic in the test species under the experimental conditions utilized.

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The ability of the DBM device test article extracts to induce chromosomal aberrations in CHO cells in the presence or absence of a rat liver homogenate metabolic activation system was evaluated per ISO 10993-3, 2003.

The DBM device test article (referenced as “Collagen membrane – EO sterile” in the report) was prepared by separately extracting 150cm² in Ham’s F-12 medium and 30 cm² in DMSO 37±1°C for 72±2 hrs. The controls used were as follows: Plain Ham’s F-12 medium and DMSO for negative controls, Mitomycin C (MMC) for the non-activated system positive control and Cyclophosphamide (CP) for the activated system positive control.

The CHO cells were treated with the Dental Barrier Membrane test article extracts in the presence and absence of metabolic activation for 3 hrs, washed with PBS and incubated with complete medium for another 21 hrs. Finally, the cells were treated with colcemid to stop metaphase for an hour before harvesting and slide preparation. The slides were scored and negative and positive responses were determined. The number and type of chromosomal aberrations was recorded.

Results: Neither of the DBM device test article extracts induced a statistically significant increase in the number of structural chromosomal aberration when compared to the negative controls in both non-activated and activated conditions. Therefore, the DBM device test article is considered non-clastogenic.

Genotoxicity Test: Rodent Bone Marrow Micronucleus Assay (Appendix 07):

The ability of the DBM device test article extracts to induce micronucleated erythrocytes in rodent bone marrow was tested per ISO 10993-3, 2003. The Dental Barrier Membrane test article (referenced as “Collagen membrane – EO sterile” in the report) was prepared by placing 180 cm² of the membrane with 30 mL of either USP 0.9% sodium chloride (NaCl) or cottonseed oil (CSO). The extraction vehicles, NaCl and CSO served as the negative controls. Cyclophosphamide (CP) served as the positive control. Ninety (90) Albino Swiss mice were used and assigned to the experimental and control groups as follows:

Table: Test Groups and Time of Sacrifice

Group	Number of Animals		Time of Sacrifice (hours)
	Extract 1: NaCl	Extract 2: CSO	
Experimental Group	10	10	24
	10	10	48
Negative Control	10	10	24
	10	10	48
Positive Control	10	N/A	24

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Results showed that the DBM device test article extracts in either NaCl or CSO did not induce a statistically significant increase in micronucleated erythrocytes as compared to the respective controls. The positive control performed as expected, thus validating the assay.

Results: Under the experimental conditions employed, the DBM device test article extracts is considered non-clastogenic.

Short Term Muscle Implantation Tests (2 weeks) (Appendix O8):

The DBM device test article (referenced as "Collagen membrane – EO sterile" in the report) was directly implanted in the paravertebral muscle tissues of New Zealand White rabbits for 2 weeks to evaluate potential to induce a local response and toxic effects. Testing was conducted based on ISO 10993-6, 2007 and implantation was done in 3 rabbits. For the control, Cartimaix Collagen membranes (another trade name for ACI-Maix) of the same size were used. Five strips of the DBM device were implanted on one side of the rabbit's paravertebral muscles and the corresponding side was implanted with the control membrane for a total of 10 implants per animal.

During the 2 week period, animals were observed daily for any clinical signs. After the 2 week period, the implants were excised for gross observations. Following formalin fixation, each of the implant sites was excised from the larger mass of tissue and examined macroscopically and scored. Histological slides of hematoxylin and eosin stained sections were prepared and examined and graded for inflammatory and healing responses. Finally, the Bioreactivity Rating is calculated based on the collected scores and indicates whether there is a difference between the test and control sites.

Results: The DBM device test article membrane implantation sites were not markedly different than the control implant sites (Bioreactivity Rating of 1.8) when implanted for 2 weeks. Therefore, the DBM device is considered non-toxic under these experimental conditions when implanted for 2 weeks.

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Short Term Muscle Implantation Tests (4 weeks) (Appendix O9):

The 4 week implantation protocol followed that of the 2 week implantation as described above. However, three (3) new rabbits were used and the implants remained in the animal for a total of 4 weeks.

Results: Results of the gross and histopathological examinations revealed that the 4 week implantation was not markedly different than the control implantation sites (Bioreactivity Rating of 2.8). Accordingly, the DBM device is deemed non-toxic under these experimental conditions when implanted for 4 weeks.

Short Term Muscle Implantation Tests (13 weeks) (Appendix O10):

The 13 week implantation protocol followed that of the 2 week implantation as described above. Three (3) rabbits (1 male and 2 female) were used and the implants remained in the animal for a total of 13 weeks.

Results: Results of the gross and histopathological examinations revealed that the 13 week implantation was not markedly different than the control implantation sites (Bioreactivity Rating of 1.0). Accordingly, the DBM device is deemed non-toxic under these experimental conditions when implanted for 13 weeks.

Short Term Muscle Implantation Tests (26 weeks) (Appendix O11):

The 26 week implantation protocol followed that of the 2 week implantation as described above. Three (3) rabbits (2 male and 1 female) were used and the implants remained in the animal for a total of 26 weeks.

Results: Results of the gross and histopathological examinations revealed that the 26 week implantation was not markedly different than the control implantation sites (Bioreactivity Rating of 0.0). Accordingly, the DBM device is deemed non-toxic under these experimental conditions when implanted for 26 weeks.

Sub-Chronic Systemic Toxicity Test (90 Day Systemic Toxicity in Rats via Intramuscular Implantation) (Appendix O12):

The DBM device (referenced in the report as “Collagen membrane – EO sterile”) was evaluated

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Manufacturing Process

An overview of the manufacturing process is provided below, including identification of the critical suppliers, specifications for reagents used, and the flow of the process.

Process Flow sheet Bulk Material Collagen (BM collagen)		
Specifications / ID	Testing activities	Process
Porcine Tissues (BM collagen) QP-SC-20035	Inspection Porcine Tissues (BM collagen) (QS-TI-1042)	BM collagen – Tissue Collection (QS-WI-0921)
Bulk Material Collagen QP-SC-10301		BM collagen – Disinfection I + II (QS-WI-0922 / QS-WI-0923)
		BM collagen – Chemical defatting (QS-WI-0924)
		BM collagen – Alkaline treatment (QS-WI-0925)
		BM collagen – Acidic treatment (QS-WI-0926)

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		BM collagen – pH adjustment (QS-WI-0927)
		BM collagen – Drying (QS-WI-0928)
	Inspection of BM collagen (QS-TI- 1043)	
		BM collagen – Batch release (QS-FO-0933)

Process Flow sheet Dental Barrier Membrane (DBM)		
Specifications / ID	Testing activities	Process
Bulk Material Collagen QP-SC-10301	Released BM collagen	
DBM – non- sterile finished product QP-SC-10302		DBM - Confectioning (QS-WI-0929 / QS-FO-0928)
		DBM - Final Inspection & Primary packaging (QS-WI-0930 / QS-FO-0928)
		DBM - Secondary packaging (QS-WI-0931 / QS-FO-0929)

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	ISO 10993-1	Biological Evaluation of Medical Devices	YES	
	ISO 10993-3	Biological evaluation of medical devices—Part 3: Tests for genotoxicity, carcinogenicity, and reproductive toxicity	YES	
	ISO 10993-4	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood - Incorporates Amendment 1: 2006	YES	
	ISO 10993-5	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity	YES	
	ISO 10993-6	Biological evaluation of medical devices—Part 6: Tests for local effects after implantation	YES	
	ISO 10993-7	Biological evaluation of medical devices — Part 7: Ethylene oxide sterilization residuals - Second Edition	YES	
	ISO 10993-10	Biological evaluation of medical devices — Part 10: Tests for irritation and skin sensitization	YES	
	ISO 10993-11	Biological evaluation of medical devices—Part 11: Tests for systemic toxicity	YES	

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	allowing bone regeneration to occur.
Application method	Quick rehydration. Space making material needed. Fixation in place not needed. Fixation with pins or sutures is possible.
Size	15x20 mm, 25x30 mm, 30x40 mm
Color	White/off white
Form	Membrane (rectangle with rounded edges)
Cross-linking	No
Resorbable	Yes
Conformability	The membrane is flexible and conforms to the defect site
Non-pyrogenic	Yes
Biocompatibility	Yes
Sterilization	Ethylene oxide
Shelf Life	3 years

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Rx only

Caution: U. S. federal law restricts this device for sale by or on the order of a physician.

Consult Instructions for Use.

Reject product if packaging damaged or opened

Information in Russian sticker

Russian sticker for final package will be added locally in Russia after approval with information below (information can be changed in accordance with Russian legislation)

Product name (in Russian)

Article:

Lot:

Manufacturer:

Place of manufacture:

Authorized representative, seller:

Sterile. Shelf life: see packaging. For single use only.

Non-pyrogenic

Non-toxic

Do not use if damaged packaging. Refer to instructions for use

ONLY FOR PROFESSIONAL USE IN STOMATOLOGY

RC # _____ date of issue _____

Label of Outer Chipboard box consist the following information:

- Order number
- Shipment number
- Barcode
- Caution about sterile content
- Cautions symbols

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Datum
13.05.2013

Certificate Number: 02444

Certificate for an animal material To be used in medical devices

Chargen No.:20035-1305

Country in which animal was slaughtered

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Geographical source of animal

Nordrhein-Westfalen

Approval Number of Slaughterhouse

ES 269

Establishment name and location

Bernhard Frenken GmbH
Paradiesstr. 19
D – 52349 Düren

Material of (state animal species)

Schwein

Age of animal

5 – 6 Monate

Nature and amount of Tissue or organ

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

Health attestation

I, the undersigned veterinarian, certify that the material described in Certificat Number 02444 was obtained from animals considered to be fit for human consumption and subject to ante- and post mortem inspection:

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Certificate Number: 1012

Certificate for an animal material
to be used in medical devices

Chargen No.: 20020-1204

Country in which animal was slaughtered
Deutschland

Bundesrepublik

Geographical source of animal

Nordrhein-Westfalen

Approval Number of Slaughterhouse

ES 1733

Establishment name and location

FVE Eschweiler eG
Indestr. 20
D-52249 Eschweiler

Material of (state animal species)

Schwein

Age of animal

5 - 6 Monate

Nature and amount of tissue or organ

2 200 g
Peritoneummembran

Packaging materials used

Kunststoff

Number of containers/packages

1

Name of veterinarian

Dr. Heyde

Collection date(s)

05.09.2012

Signed _____
(Responsible person)

Health attestation

I, the undersigned veterinarian, certify that the material described in Certificate Number 1012 was obtained from animals considered to be fit for human consumption and subject to ante- and post-mortem inspection:

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herewith confirms that the company

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has introduced, applies and maintains a Quality Management System in the area of:

Sterilization of medical devices with Ethylenoxide considering the
DIN EN ISO 11135:2014 standard including corresponding
microbiological and chemical inprocess and final controls


The compliance of the Quality Management System with the requirements of the below mentioned standard was verified by an audit:

EN ISO 13485:2012 + AC:2012

This certificate is valid until: 31 March 2019

Report No.: 2100FS17F
Process No.: QS – 2100
Certificate No.: 2100GB438160408B

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2. REVISION RECORD

Rev. A	First release.
Rev. B	The need to update the first revision resulted from the following: • The decision to put Remaix – ortho development on hold (QP-MT-03-23Jun08) • The design of the product and the manufacturing process have been finalized (QP-PF-03001-C) • The outcome of a Design Process audit held by qcmed on 04.Dec.08 as part of the internal auditing program. As presented in Audit-report-Mat-081204, 2 recommendations for improvement were made. These were incorporated in this Risk Analysis. • The outcome of the clinical evaluation of the Remaix membrane (QP-CR-03001-A).
Rev. C	Product nomenclature revised according to Engineering Change Notice QP-EC-12002. Qualifications of personnel involved updated.
Rev. D	Document release section improved to prevent that originators check and release their own work. creos TM xeno.protect added to the product family (QP-EC-13003). Formal update of the FMECA to include the ALAP instead of ALARP definition according to EN ISO 14971:2012.
Rev. E	Inclusion of results of the DBM shelf life validation study data. Inclusion of the risk assessment related to the introduction of the rotary die cutter in section 12.5 (see QP-RA-03006-A). Update of section 12.2 and 12.3 since the control measures implemented for risk 12.5 also improve the control of risks 12.2 and 12.3. References to clinical evaluations updated
Rev. F	Formal update of the table in section 6.3 since the EO dot is no longer applied on the second sterile barrier system since the introduction of bulk sterilization.

3. IDENTIFICATION AND QUALIFICATION OF PERSONNEL INVOLVED

This risk analysis was performed by the Matricel representatives indicated in the table below. In addition this table indicates the qualification and expertise of the people involved in this risk analysis.

Name	Qualification	Expertise*
Ingo Heschel	Mechanical Engineer by training, previous Head of R&D Matricel, >23 years experience in medical device research and development, clinical expert ACI-Maix and DBM.	R, E, F, C.
Leon Olde Damink	Chemical Engineer by training, Head of R&D, Head of RA-QM, previous Head of Manufacturing Matricel, >18 years experience in medical device development and manufacturing, clinical expert ACI-Maix and DBM.	R, M, F, C.
Lizette Eummelen	Chemical Engineer by training, Head of QC, >13 years experience in medical device development and manufacturing.	M.
Frank Schügner	Mechanical Engineer by training, Head of Engineering, >13 years experience in engineering.	E, M.
Hans Leemhuis	Mechanical Engineer by training, Head of Manufacturing, >5 years experience in engineering.	E,M.

*R=Risk management expert, E=Engineering/Construction, M=Manufacturing, F=Function of device, C=Clinical use of device.

4. REFERENCES

QP-DC-03001: Design Concept Remaix membrane for bone graft containment
QP-RP-03001: Risk Management Plan Remaix membrane for bone graft containment.
QP-PF-03001: Process Flow sheet Bulk Material Collagen manufacturing and testing
QP-PF-03002: Process Flow sheet Dental Barrier Membrane manufacturing and testing

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

	Risk Analysis	Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971	Release: 07.Mar.16
		Page: 3 of 22

QP-RA-03001: Dental Barrier Membrane - Risk Analysis EN ISO 22442-1
 QP-RA-03006: Confectioning of membranes with rotary die cutter
 QP-CR-03001: Remail membrane - Evaluation of clinical data
 QP-CR-03002: Matricel Dental Barrier Membrane – CE-certified product - Evaluation of clinical data
 QP-CR-03003: Matricel Dental Barrier Membrane – Health Canada submission - Evaluation of clinical data
 EN ISO 14971: Medical devices - Application of risk management to medical devices.
 EN ISO 22442-1: Medical devices utilizing animal tissues and their derivatives - Part 1: Application of risk management.

5. DEFINITIONS

The following definitions apply to the risk analysis described in this document:

Term	Definition
FMECA	Failure mode, effects and criticality analysis.
Risk analysis	Systematic use of available information to identify hazards and to estimate the risk.
Harm	Physical injury or damage to the health of people, or damage to property or the environment.
Hazard	Potential cause of harm.
Hazardous situation	Circumstance in which people, property, or the environment are exposed to one or more hazard(s).
Failure mode	(Sequence of) event(s) that leads to a hazard or hazardous situation.
Risk	Combination of the probability of occurrence of harm and the severity of that harm.
Numerical Risk Value	Value calculated from the FMECA that allows an assessment of risk acceptability.
Safety	Freedom from unacceptable risk.
Risk control	Corrective actions needed to bring a risk to an acceptable level.
Verification	Confirmation, through the provision of objective evidence, that specified requirements have been fulfilled.

6. USE, DESCRIPTION AND MANUFACTURING PROCESS OF THE PRODUCT

6.1 Description of the product, indications.

A detailed description of the product and its use, users and use process can be found in the Design Concept document (QP-DC-03001) and the Risk Management Plan (QP-RP-03001) for this project.

6.2 Description of the production process

The DBM product is prepared from Bulk Material Collagen (BM collagen), the purified collagen used by Matricel to manufacture collagen products. Please refer to document QP-PF-03001 for a detailed description of the BM collagen production process. Please refer to document QP-PF-03002 for a detailed description of the DBM product production process.

6.3 Functional description of the components

All components of the DBM membrane and packaging referred to in the risk analysis matrix (see Attachment 2) are listed here with a comprehensive description of their function.

	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 4 of 22

Component	Function
DBM membrane	Collagen membrane for the containment of bone graft material.
First sterile barrier system	Provides a barrier for microorganisms to keep the DBM membrane sterile inside. Allows the sterilization of the device.
Second sterile barrier system	Holds the inner packaging. Has a sterile inner surface. Carries the product label. Allows the sterilization of the device.
Product label	Provides product- and manufacturer information.
Product box	Contains the packaged product and Instructions for Use. Provides product-handling and manufacturer information.
Product box label	Provides product information.
Instructions for Use	Provides information to the user.
Documentation label	Provides product- and manufacturer information.
Shipping box	Contains the packaged product and the EO dot during shipment to the sterilizer.
EO dot	Indicates that EO sterilization has been performed.
Customer shipping box	Contains the product during shipment to the customer.

7 RISK ANALYSIS

The failure mode and criticality analysis (FMECA) as described in QS-SP-0403 is used as a tool employed for identifying and evaluating hazards and risks associated with a product and/or process. Compared to the guidance presented in the SOP, the quantitative assessment for Occurrence and Detection has been adapted to the expected annual sales volume of around 10.000 products. This will allow a meaningful comparison with future manufacturing trend analysis data and Post Market Surveillance data.

The results of the FMECA are presented in the Risk Analysis Matrix provided as Attachment 2 to this report. This FMECA covers the risk related to BM collagen as well as the DBM product.

7.1 Guidance on the contents of the Risk Analysis Matrix

Column 1. Identification of characteristics that could impact safety. EN ISO 14971 requires that the manufacturer identifies those characteristics of the medical device that could affect safety. The guidance provided in Annex C is used as an aid in identifying all the characteristics of the medical device that could affect safety. These questions were answered from the point of view of all the individuals involved (e.g., users, maintainers, patients, etc.). The list is not exhaustive, or representative of all medical devices, and the reader is cautioned to add questions that can have applicability to the particular medical device and to skip questions that are not relevant to the particular medical device (rationale provided in Attachment 2). The reader is also cautioned to not only consider each question on its own but also in relation to each other.

Column 2: Running number for identifying each hazard or hazardous situation.

Column 3: Hazard or hazardous situation as a result of a failure mode.

Column 4: A description of the failure mode that must occur to cause the hazard or hazardous situation. A failure mode can be the result of a single event or a sequence of events.

Column 5: The number of events described in column 4 necessary for the failure mode to occur is indicated (n=0,1,2,...). The highest risk is correlated to failure modes that could occur without any event (i.e. n=0), in other words occurrence of a hazard or hazardous situation although the product is „perfect“.

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 5 of 22

Column 6: Severity: Indicator for the severity of the potential harm resulting from a hazard or hazardous situation to the users, maintainers, patients, etc. Score is 1 (one, best) to 10 (ten, worst). Detailed definitions are provided in Attachment 1.

Column 7: Occurrence: Likelihood that a failure mode will occur. Score is 1 (one, best) to 10 (ten, worst). Detailed definitions are provided in Attachment 1.

Column 8: Detection: Likelihood that a failure mode will be identified in case it occurs. Score is 1 (one, best) to 10 (ten, worst). Detailed definitions are provided in Attachment 1.

Column 9: Numerical Risk Value (NRV): Multiplication of score in columns 6, 7 and 8 to obtain an overall numerical value for the risk. Score is 1 (one, best) to 1000 (thousand, worst). The colour of the numerical risk value is indicating the acceptance of the risk as defined under 7.2.

Column 10: Risk Control/Verification. If NRV can not be calculated (for instance early in the design process), column to indicate the control measures and verification activities needed to allow calculation. If NRV is unacceptable, column to indicate the risk control measures (corrective actions) needed to achieve an acceptable NRV (note that the implemented risk control measures shall be reviewed with regard to the introduction of new hazards or hazardous situations as described in the Risk Management Plan). If NRV is acceptable, column to indicate the verification activities performed or controls put in place.

Note that in case the NRV is not acceptable and corrective actions are needed to achieve an acceptable NRV, the risk analysis matrix provided in attachment 2 will be expanded with additional columns to document the verification of the effectiveness of the corrective action and to evaluate the residual risk after implementation of the corrective action.

7.2 Acceptance criteria for risks

Each risk is evaluated according to the three parameters: occurrence, severity and detection. In general, there are three possible ranges for the numerical risk value. However, the final acceptance of each risk is the ultimate decision of the responsible manager. For risks that are completely under control and responsibility of Matricel, a value in Range 1 is acceptable. For risks that are out of the control of Matricel, such as e.g. storage and handling at the customer, a value in Range 2 is acceptable. For all Risks falling into Range 3, risk control measures need to be implemented followed by a renewed assessment of the risks resulting from the implemented measures.

It is recognized, that the risk associated with a device, its components, production and use, can not fully be evaluated by presenting only a numerical value. So extensive engineering and medical knowledge are used to evaluate the acceptance for each individual risk.

Range	Risk acceptance	Numerical Risk Value	Justification
Range 1	generally accepted	1 to 100	NRV values of maximally 100 will guarantee that: A. Failure modes with high Occurrence (Occurrence range 8-10) and severe effects (Severity range 7-10) will only be accepted if they have a very high likelihood of detection (Likelihood of Detection 1) before they occur. B. Failure modes with severe effects (Severity range 7-10) and a low likelihood of detection (Likelihood of Detection 8-10) will only be accepted if they have a very low Occurrence (Occurrence 1).
Range 2	as low as possible	101 to 150	Failure modes resulting in NRV-values in this range will only be accepted if they result in hazards with a light to moderate severity.
Range 3	generally not accepted	151 to 1000	Not applicable.

	Risk Analysis	Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971	Release: 07.Mar.16
		Page: 6 of 22

7.3 Risk Analysis Matrix

The risk analysis matrix is presented as Attachment 2 to this report. The outcome of the risk analysis as presented in column 'Risk' is also represented in the colours indicated above. For column 'Risk control/Verification', in case corrective actions are needed for risk control, they are identified in red, in case verification activities are required they are identified in light blue. As a consequence, at the end of a successful risk management process all NRV values should be coloured green.

8 STATEMENT ON RISK ACCEPTANCE

All reasonably foreseeable risks associated with the design, manufacture and use of the Dental Barrier Membrane have been taken into account for this risk analysis. The evaluation has been performed to the company's best knowledge.

The outcome of this risk analysis is that the risks associated with all hazardous situations taken into consideration are controlled to an acceptable level.

Verification and validation activities (performed voluntarily or as required by European and international standards) that have been performed to support the risk assessment for each hazardous situation, have been closed successfully.

The clinical benefit of the DBM product has been assessed in the current versions of the clinical evaluation documents.

Weighing the clinical benefits presented in the current versions of the clinical evaluation documents and the risks assessed in this risk analysis, it is concluded that the benefits of using the DBM product in dental bone augmentation procedures outweigh the risks associated with its use.

This benefit/risk assessment is made to the best of our knowledge.

	Risk Analysis	Document: QP-RA-03002-F
		Release: 07.Mar.16
	Dental Barrier Membrane Risk Analysis EN ISO 14971	Page: 7 of 22

ATTACHMENT 1: Overview of classifications used for risk analysis.

Column 6: Severity.

Denominator	Description of harm	Factor
Light harm potential	Transient adverse effect on body functions	1
	Light injury	2
	Light impairment of the health status	3
Moderate harm potential	Short-term impairment of body functions (days)	4
	Moderate injury	5
	Moderate deterioration in health status	6
	Permanent impairment of body functions	7
High harm potential	Severe deterioration of health status	8
	Life-threatening condition	9
	Occurrence of the failure will most probably be followed by death	10

Column 7: Occurrence.

Denominator	Qualitative assessment	Quantitative assessment	Factor
Very low	Event is not probable.	< 1 in 10.000	1
Remote	Event may occur seldom.	1 in 10.000	2
		1 in 2000	3
Occasional	Event has been found in few cases.	1 in 500	4
		1 in 200	5
		1 in 100	6
Probable	Event will probably occur in some cases.	1 in 50	7
		1 in 20	8
Expected	Event can almost not be avoided.	1 in 10	9
		1 in 2	10

Column 8: Detection.

Denominator	Qualitative assessment	Quantitative assessment	Factor
Very high	Process controls will surely reveal a potential event.	< 1 in 10.000	1
		1 in 10.000	2
High	Process controls have a good chance to reveal a potential event.	1 in 2000	3
		1 in 500	4
Moderate	Process controls will probably reveal a potential event.	1 in 200	5
		1 in 100	6
Low	Process controls may not reveal a potential event.	1 in 50	7
		1 in 20	8
Very low	Process controls can hardly reveal a potential event.	1 in 10	9
Improbable	Process controls can not reveal a potential event.	1 in 2	10

	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 8 of 22

ATTACHMENT 2: Risk Analysis Matrix for the Dental Barrier Membrane for dental bone augmentation.

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
1. What is the intended use and how is the medical device to be used? DBM product provides barrier to protect bone grafted defect from ingrowing connective tissue. Does not provide life support. Can be easily removed from the implant site in case of complications.	1.1	Failure of treatment due to lack of initial barrier function.	Impaired integrity of DBM product due to puncture or rupture of the starting materials during harvesting.	1	6	2	2		100 % verification of freedom of defects during product inspection.
	1.2	Failure of treatment due to lack of initial barrier function.	Impaired integrity of DBM product due to puncture or rupture of the product during manufacturing.	1	6	2	2		100 % verification of freedom of defects during product inspection.
	1.3	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to impaired collagen nativity resulting from incorrect storage after harvesting.	1	6	2	2		Maximum storage time after harvesting validated. In-process monitoring in place.
	1.4	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to impaired collagen nativity resulting from incorrect chemical concentrations during manufacturing.	1	6	2	2		In-process monitoring in place. Manufacturing process validated.
	1.5	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to impaired collagen nativity resulting from incorrect processing time during manufacturing.	1	6	2	2		In-process monitoring in place. Manufacturing process validated.
	1.6	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to impaired collagen nativity resulting from excessive temperatures during manufacturing.	1	6	2	2		In-process monitoring in place. Manufacturing process validated.
	1.7	Failure of treatment due to early loss of barrier function after implantation.	Too quick biodegradation due to insufficient thickness of the product.	1	6	2	2		100 % verification of product weight during product inspection.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 9 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
	1.8	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to infection of implant site.	0	6	6	2		Warning in IFU not to implant in areas with active infection.
	1.9	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to early exposure as a result of soft tissue dehiscence.	0	6	8	2		Soft tissue dehiscence is a common, minor complication. Can be due to surgical as well as membrane related factors. Based on bio-equivalence data, DBM product is expected to have similar complication rates as the current market leading product. Surgical related factors (i.e. tension free closure) are included in the IFU.
2. Is the medical device intended to be implanted? Procedure will be performed with patients of all ages, weight, and physical activity. DBM product is implanted in contact with jaw bone in	2.1	Incompatibility with host tissue resulting in unwanted tissue reaction.	Ineffective purification resulting in the presence of non-collagenous impurities in the product that could potentially include immunogens.	1	8	2	2		In-process monitoring in place. Manufacturing process validated.
	2.2	Incompatibility with host tissue resulting in unwanted tissue reaction.	Ineffective removal of manufacturing remnants resulting in the presence of excessive contamination.	1	4	2	2		In-process monitoring in place. Manufacturing process validated.
	2.3	Incompatibility with host tissue resulting in unwanted tissue reaction.	Impaired collagen nativity resulting from incorrect manufacturing conditions.	1	3	2	2		In-process monitoring in place. Manufacturing process validated.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 10 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
patients with normal health requiring local bone regeneration.	2.4	Incompatibility with host tissue resulting in unwanted tissue reaction.	Long term contact with host tissue (more than 30 days) resulting in unwanted tissue reaction.	0	3	2	10		Biocompatibility according to ISO 10993 series assessed.
	2.5	Incompatibility with host tissue resulting in unwanted tissue reaction.	Single exposure to product resulting in unwanted tissue reaction related to presence of collagen-related immunogens.	0	4	2	10		Biocompatibility according to ISO 10993 series assessed. Warning in Instructions for Use.
	2.6	Incompatibility with host tissue resulting in unwanted tissue reaction.	Multiple exposure to product resulting in unwanted tissue reaction related to presence of collagen-related immunogens.	0	4	1	10		Biocompatibility according to ISO 10993 series assessed. Warning in Instructions for Use.
3. What materials or components are incorporated in the medical device or are used with, or are in contact with, the medical device? Are there unwanted outputs of energy or substances?	3.1	Unwanted biological reaction due to presence of parasites.	Ineffective control of health status of animal followed by survival of parasites during manufacturing.	2	9	1	1		Harvesting compliant to EN ISO 22442-2. Disinfection of fresh membranes validated. Manufacturing process validated. Sterilization process validated. See current version of QP-RA-03001 for extensive assessment.
DBM product is composed of purified collagen	3.2	Unwanted biological reaction due to presence of unclassified pathogenic entities.	Survival of unclassified pathogenic entities during manufacturing.	1	9	1	10		Active surveillance of development of zoonotic treats. See current version of QP-RA-03001 for extensive assessment.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 11 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
and elastin. Remnants of animal sourcing or the manufacturing process may be present. DBM product does not emit energy.	3.3	Unwanted biological reaction due to presence of excessive bioburden content.	Ineffective control of bioburden contamination, resulting in too high bioburden content, followed by routine sterilization, resulting in non-sterile product.	1	9	2	1		Disinfection of fresh membranes validated. Bioburden level is under control. Sterilization validated. See current version of QP-RA-03001 for extensive assessment.
	3.4	Unwanted biological reaction due to presence of excessive bioburden content.	Ineffective control of bioburden contamination, resulting in too high bioburden content, followed by ineffective sterilization, resulting in non-sterile product.	2	9	1	1		Disinfection of fresh membranes validated. Bioburden level is under control. Sterilization validation. See current version of QP-RA-03001 for extensive assessment.
	3.5	Unwanted biological reaction due to presence of viral contamination.	Ineffective control of health status of animal and survival of viral contamination during manufacturing.	2	10	1	10		Harvesting compliant to EN ISO 22442-2. Efficiency of manufacturing process for viral inactivation validated. See current version of QP-RA-03001 for extensive assessment.
	3.6	Unwanted biological reaction due to presence of TSE transmission.	Ineffective control of health status of animal and survival of TSE contamination during manufacturing (<u>hypothetical case</u>).	x	10	1	10		Harvesting compliant to EN ISO 22442-2. Prior inactivation step incorporated in process. See current version of QP-RA-03001 for extensive assessment.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 12 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
	3.7	Unwanted biological reaction due to presence of pyrogens.	Ineffective control of bioburden during manufacturing followed by effective sterilization.	1	9	3	3	81	Processing water is monitored. Process validation. Endotoxin testing for every batch. See current version of QP-RA-03001 for extensive assessment.
	3.8	Unwanted biological reaction due to presence of immunogens.	Ineffective purification resulting in the presence of non-collagenous impurities in the product that could potentially include immunogens.	1	8	2	2	32	In-process monitoring in place. Manufacturing process validated.
	3.9	Unwanted biological reaction due to remnants of the manufacturing process.	Ineffective removal of acetone used during the manufacturing process.	1	4	2	1	3	Pharmaceutical grade acetone used if remnants can end in product. In-process monitoring in place. Acetone residues tested for every batch. Manufacturing process validated.
	3.10	Unwanted biological reaction due to remnants of the manufacturing process.	Ineffective removal of residual salts generated during the manufacturing process.	1	3	2	1	3	Pharmaceutical grade chemicals used if remnants can end in product. In-process monitoring in place. Ash residue tested for every batch. Manufacturing process validated.
	3.11	Unwanted biological reaction due to remnants of the manufacturing process.	Excessive acidity of the product.	1	6	2	1	12	In-process monitoring in place. pH tested for every batch. Manufacturing process validated.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 13 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
	3.12	Unwanted biological reaction due to substances extracted from equipment used during the manufacturing process.	Extraction of compounds from equipment and auxiliary materials during contact with product and production fluids during manufacturing process.	1	8	3	3		Minimum quality of materials with product contact is "Food grade". Risk analysis for genotoxicity performed. Genotoxicity testing performed.
	3.13	Unwanted biological reaction due to degradation products of the membrane.	Ineffective control of manufacturing conditions resulting in (partial) denaturation of collagen.	n	3	2	2		In-process monitoring in place. Denaturation temperature measured for every batch. Manufacturing process validated.
	3.14	Unwanted biological reactions due to remnants of the cleaning agents used.	Ineffective control of removal of cleaning agents.	1	4	2	2		Cleaning validation performed.
	3.15	Unwanted biological reactions due to foreign particulate present on the product.	Ineffective control of environmental conditions during manufacturing. Failure to identify presence of foreign particulate during product inspection.	2	8	3	2		Yearly clean room qualification. Frequent clean room monitoring (viable and non-viable particles). 100 % verification of freedom of foreign particulate during product inspection. Compatibility of equipment with clean room environment demonstrated during DQ-IQ.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 14 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
4. Is the medical device supplied sterile or intended to be sterilized by the user, or are other microbiological controls applicable? DBM product is delivered sterile and for single use only.	4.1	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to too high bioburden before sterilization.	1	9	1	1	1	Process monitoring. Process validation. Bioburden monitoring for every batch. Shelf life of non-sterile product validated.
	4.2	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to ineffective sterilisation.	1	9	1	1	1	Sterilization validated. Product release for every batch based on parametric as well as bio-indicator release.
	4.3	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to expiration of shelf life.	1	9	1	1	1	Warning on product label and Instruction for Use.
	4.4	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to damaged inner packaging + damaged outer packaging.	2	9	1	2	18	Warning on product label and Instruction for Use.
	4.5	Unwanted biological reaction due to un-sterile product.	Deterioration of properties due to re-sterilisation by the end user.	0	9	1	10	30	Warning on product label and Instruction for Use.
	4.6	Unwanted biological reaction due to presence of pyrogens.	Ineffective control of bioburden contamination followed by effective sterilization.	1	9	3	3	31	Processing water is monitored. Process validation. Endotoxin testing for every batch. See current version of QP-RA-03001 for extensive assessment.
5. Is the medical device susceptible to environmental influences?	5.1	Failure of treatment due to lack of initial barrier function.	Impaired integrity of DBM product due to damage of the product during transport through the first and second sterile barrier system and the product box (puncture or rupture).	3	6	3	3	44	Shipping validation data available confirm protective properties of packaging system. Warning on product label.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 15 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
DBM product is susceptible to excessive mechanical force, high temperature in combination with high humidity and excessive exposure to light. The protective package is a Tyvek lid that is permeable for fluids.	5.2	Failure of treatment due to early loss of barrier function after implantation.	Increased degradation due to impaired collagen nativity resulting from incorrect storage after harvesting. High temperature in combination with high humidity. Ineffective control of shipping conditions resulting in (partial) denaturation of collagen.	n	6	2	5	50	Storage conditions after harvesting specified and validated. Storage conditions of end product specified in Instruction for Use and validated during shelf life study.
	5.3	Failure of treatment due to early loss of barrier function after implantation.	Increased degradation due to impaired collagen nativity resulting from incorrect storage after harvesting. Exposure to light. Ineffective control of shipping conditions resulting in (partial) denaturation of collagen.	n	6	1	5	30	Product packaged in product box. Recommended storage conditions specified in Instructions for Use.
	5.4	Failure of treatment due to impaired handling properties.	Membrane has been exposed to liquid during distribution. Subsequent drying results in inflexible product.	1	6	2	3	10	Recommended storage conditions specified and warning in Instructions for Use.
6. Does the medical device have a restricted shelf-life?	6.1	Failure of treatment.	Impaired properties due to expiration of shelf life	1	6	1	1	10	Warning not to use expired product in Instructions for Use.
DBM product has a limited shelf life. There are no risks related to the disposal of remnants of the product when not	6.2	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to expiration of shelf life.	1	9	1	1	10	Warning not to use expired product in Instructions for Use.
	6.3	Misuse of product resulting in failure of treatment.	Impaired readability of the information provided resulting from events during shipment.	1	6	2	2	10	Shelf life validation data available. Shipping validation data available confirm stable properties of packaging system.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 16 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
used.	6.4	Misuse of product resulting in failure of treatment.	Loss of label during shelf life of the product.	1	6	2	2		Shelf life validation data available. Shipping validation data available confirm stable properties of packaging system.
7. To what mechanical forces will the medical device be subjected?	7.1	Failure of treatment due to lack of initial barrier function.	Impaired integrity of DBM product due to damage of the product during clinical application (puncture or rupture).	1	6	2	3		Mechanical properties of product verified to be compliant with design input requirements.
DBM product will be applied with standard surgical tools like tweezers. Typically DBM product will not be fixed but in case needed, DBM product may be sutured. DBM product will be exposed to typical transport related forces during distribution.	7.2	Failure of treatment due to lack of initial barrier function.	Excessive forces during storage and distribution.	1	6	3	3		Shelf life validation data available. Shipping validation data available confirm stable properties of packaging system.
	7.3	Failure of treatment due to ineffective fixation.	Mechanical properties of DBM product insufficient to allow preferred fixation technique.	1	6	2	3		Mechanical properties of product verified to be compliant with design input requirements.
	7.4	Failure of treatment due to ineffective fixation.	Mechanical properties of DBM product insufficient to allow preferred fixation technique due to damage of the product during distribution.	1	6	3	3		Shelf life validation data available. Shipping validation data available confirm stable properties of packaging system.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 17 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
8. What determines the lifetime of the medical device? The shelf life of DBM product is determined by the collagen stability and the packaging stability.	8.1	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to impaired collagen nativity due to expiration of shelf life.	1	6	1	1	1	Shelf life validation data available.
	8.2	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to impaired packaging properties due to expiration of shelf life.	1	9	3	3	31	Shelf life validation data available. Shipping validation data available confirm stable properties of packaging system.
9. Is the medical device intended for single use? DBM product is for single use only and should not be re-sterilized.	9.1	Failure of treatment due to early loss of barrier function after implantation.	Increased degradation due to impaired collagen nativity resulting from re-sterilisation by the end user.	0	6	1	10	30	Warning on product label and Instruction for Use.
	9.2	Unwanted biological reaction due to un-sterile product.	Re-use of remnants of the product after opening of the package for the next patient.	0	9	1	10	30	Warning on product label and Instruction for Use.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 18 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
10. Does installation or use of the medical device require special training or special skills? The use of barrier membranes like DBM product is well established technique.	10.1	Failure of treatment due to inappropriate use of the product.	Ineffectiveness of training.	1	6	2	2	32	The use of barrier membranes in dental bone augmentation is a very well established technique. Information provided on product label and in Instruction for Use.
11. How will information for safe use be provided? DBM product will be packaged individually. Every package will include Instruction for Use and sold directly to end user.	11.1	Failure of treatment due to inappropriate use of the product.	Missing Instruction for Use.	1	6	2	2	24	In-process control and Finished product inspection in place.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 19 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
12. Will new manufacturing processes need to be established or introduced? The DBM product manufacturing process could include an optimized drying method, an optimized trimming method, an automated cutting method and the use of a new packaging configuration.	12.1	Unwanted biological reaction due to substances extracted from equipment used during the manufacturing process.	Optimized drying method. Extraction of compounds from new equipment during contact with product and production fluids during manufacturing process.	1	8	3	3		Minimum quality of materials with product contact is "Food grade". Compatibility of equipment with product demonstrated during DQ-IQ.
	12.2	Unwanted biological reactions due to foreign particulate present on the product.	Generation of foreign particles as a result of the optimized trimming method.	1	8	3	1		Compatibility of equipment with clean room environment and product demonstrated during DQ-IQ. 100 % verification of freedom of foreign particulate during product inspection. Sampling for heavy metals testing after trimming of product.
	12.3	Unwanted biological reactions due to foreign particulate present on the product.	Generation of foreign particles as a result of the automated cutting method.	1	8	3	1		Compatibility of equipment with clean room environment and product demonstrated during DQ-IQ-OQ. 100 % verification of freedom of foreign particulate during product inspection. Sampling for heavy metals testing after die-cutting of product.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 20 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
	12.4	Unwanted biological reaction due to un-sterile product.	Product un-sterile due to impaired packaging properties due to new packaging configuration.	1	9	3	3	■	Qualification of packaging system performed. Shipping validation data available confirm stable properties of packaging system.
	12.5	Unwanted biological reactions due to foreign particulate present on the product.	Generation of foreign particles as a result of the rotary die cutting method.	1	8	3	1	■	Compatibility of equipment with clean room environment and product demonstrated during DQ-IQ. Warnings in equipment instruction. 100 % verification of freedom of foreign particulate during product inspection. Sampling for heavy metals testing after die-cutting of product.
13. In what way(s) might the medical device be deliberately misused?	13.1	Failure of treatment due to early loss of barrier function after implantation.	Impairment of nativity of the collagen. Neglect of manufacturer's instruction related to re-sterilisation.	0	6	1	10	■	Warning on product label and Instructions for Use.
	13.2	Unwanted biological reaction due to un-sterile product.	Impairment of sterility of the product Neglect of manufacturers instruction related to single use only.	0	9	1	10	■	Warning on product label and Instructions for Use.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis		Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971		Release: 07.Mar.16
			Page: 21 of 22

1	2	3	4	5	6	7	8	9	10
Identification of characteristics that could impact safety	#	Hazardous situation	Failure mode	n	Sev.	Occ.	Det.	NRV	Risk Control/ Verification
14. Does the use of the medical device depend on essential performance requirements?	14.1	Failure of treatment due to lack of initial barrier function.	Impaired integrity of DBM product due to puncture or rupture of the membrane.	n	6	2	2	14	All controls in place to prevent the distribution of DBM product with an impaired integrity.
	14.2	Failure of treatment due to early loss of barrier function after implantation.	Increased biodegradation due to impaired collagen nativity resulting from incorrect storage after harvesting.	n	6	2	2	14	All controls in place to prevent the distribution of DBM product with an impaired nativity.
	14.3	Failure of treatment due to early loss of barrier function after implantation.	Persistence of membrane too short due to insufficient thickness of the product.	1	6	2	2	14	All controls in place to prevent the distribution of DBM product with insufficient thickness.
	14.4	Failure of treatment due to lack of initial barrier function.	Impaired integrity of DBM product due to damage of the product during clinical application (puncture or rupture).	1	6	2	3	30	Mechanical properties of product verified to be compliant with design input requirements.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

	Risk Analysis	Document: QP-RA-03002-F
	Dental Barrier Membrane Risk Analysis EN ISO 14971	Release: 07.Mar.16
		Page: 22 of 22

ATTACHMENT 3: Overview of questions listed in Annex C of EN ISO 14971 that are not applicable to Dental Barrier Membrane and rationale

#	Questions not applicable	Rationale
1	Is energy delivered to or extracted from the patient?	DBM product is a non-active implantable medical device.
2	Are biological materials processed by the medical device for subsequent reuse, transfusion or transplantation?	DBM product does not process biological materials.
3	Is the medical device intended to be sterilized by the user, or are other microbiological controls applicable?	DBM product is not intended to be sterilized or re-sterilized by the user.
4	Is the medical device intended to be routinely cleaned and disinfected by the user?	DBM product is not intended to be routinely cleaned and disinfected by the user.
5	Is the medical device intended to modify the patient environment?	DBM product is not intended to modify the environment with respect to temperature, humidity, atmospheric gas composition, pressure, and light.
6	Does the medical device influence the environment?	DBM product is a non-active implantable medical device.
7	Are measurements taken?	DBM product does not have a measuring function.
8	Is the medical device interpretative?	DBM product does not have a measuring function.
9	Is maintenance or calibration necessary?	DBM product does not have a measuring function.
10	Is the medical device intended for use in conjunction with medicines or other medical technologies? Are there essential consumables or accessories associated with the medical device?	DBM product is not intended to be used in conjunction with medicines or other essential medical technologies/accessories.
11	Does the medical device contain software?	DBM product is a non-active implantable medical device.
12	Is safe decommissioning or disposal of the medical device necessary?	DBM product and the packaging materials used do not contain toxic or hazardous material and can be disposed of as normal garbage.
13	Is successful application of the medical device critically dependent on human factors such as the user interface?	DBM product is a non-active implantable medical device.
14	Does the medical device use an alarm system?	DBM product is a non-active implantable medical device.
15	Does the medical device hold data critical to patient care?	DBM product is a non-active implantable medical device.
16	Is the medical device intended to be mobile or portable?	DBM product is a non-active implantable medical device.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03002\QP-RA-03002-F.docx

matricel	Risk Analysis	Document: QP-RA-03005-B
		Release: 28.Feb.17
	Viral safety of Dental Barrier Membrane	Page: 1 of 20

PREPARED BY (QM): A. Bühren DATE: 28 Feb 17

PREPARED BY (Dr Gregersen): Gregersen DATE: 28. Febr. 17

CHECKED BY (SO): F. Schüßner DATE: 28 Feb 17

APPROVED BY (RA-QM): H. Olde Damink DATE: 28 Feb 17

1 SCOPE

This risk analysis is part of the risk management activities for the Matricel Dental Barrier Membrane (DBM) product family. This report provides an overview of the data available to demonstrate the viral safety of the DBM product as required by ISO 22442-3.

Note: This is revision B of this document. In this document the following content was updated, in close cooperation with Jens-Peter Gregersen DVM, PhD, Virologist, and Biomedicines Consultant:

- Section 3: Updated with all new references.
- Section 4: New section to describe literature review methodology, provided by Dr Gregersen.
- Section 5: Complete rework, provided by Dr Gregersen.
- Section 6: Complete rework, provided by Dr Gregersen.
- Section 10: Addition of new experimental data for viral clearance of PPV during acidic treatment, provided by Matricel.
- Section 11: Addition of new experimental data for viral clearance of MuLV, PRV, PPV, and BVDV during ethylene oxide gas sterilization, provided by Matricel.
- Section 12: Update of the discussion and critical evaluation of the overall residual risk acceptability based on the complete rework of sections 4 and 5 and the new experimental data, by Matricel and Dr Gregersen.

2 PURPOSE

The risks related to potential viral contamination is an important aspect in the assessment of the safety of the DBM product family. This document evaluates the potential viral contaminants of the porcine raw material by virus families and virus types and provides an estimate of the likely viral contamination level that might be expected under worst-case conditions. In addition, this document provides a summary of the experimental viral clearance data available for different steps of the manufacturing process of Bulk Material Collagen (BM collagen) as the starting material for DBM manufacturing, and the final sterilization process of DBM. This document concludes with a discussion and critical evaluation of the overall residual risk acceptability by comparing potential viral contamination levels with the experimental data on virus reduction.

The following sub-processes in the BM collagen and DBM manufacturing process have been experimentally validated for viral inactivation effectiveness in order to assess the viral safety risks.

	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 2 of 20

Process	Step	Experimental data	Section
BM collagen	Chemical defatting	NewLab study (PRV, PPV)	8
		Lancaster Labs study (XMuLV, PRV, Reo-3, PPV, BVDV)	9
BM collagen	Alkaline treatment	NewLab study (PRV, PPV)	8
		Lancaster Labs study (XMuLV, PRV, Reo-3, PPV, BVDV)	9
BM collagen	Acidic treatment	Charles River Labs study 2 (PPV) (added, new study data)	10
DBM	EO sterilization	Charles River Labs study 1 (MuLV, PRV, PPV, BVDV) (added, new study data)	11

3 REFERENCES

- ISO 22442-1:2015 Medical devices utilizing animal tissues and their derivatives — Part 1: Application of risk management.
- ISO 22442-3:2007 Medical devices utilizing animal tissues and their derivatives — Part 3: Validation of the elimination and/or inactivation of viruses and transmissible spongiform encephalopathy (TSE) agents
- QP-PF-03001: Manufacturing process, in-process monitoring and product testing for Bulk Material Collagen.
- QP-PF-03002: Manufacturing process, in-process monitoring and product testing for Dental Barrier Membrane.
- QP-VR-03009: Matricel Collagen Membrane - Worst case EO treatment during viral clearance study
- Study # 2/05.104: Validation of the inactivation of Pseudorabies Virus, and Porcine Parvovirus by treatment with Acetone, NaOH, and gamma irradiation in the manufacturing process of collagen under worst case conditions. Newlab BioQuality AG, 2006.
- Study # NG-1292703: Virus clearance study for porcine peritoneal tissue for the chemical defatting and alkaline treatment process steps - Xenotropic Murine Leukemia Virus (XMuLV). Eurofins Lancaster Labs, 2012.
- Study # NG-1292704: Virus clearance study for porcine peritoneal tissue for the chemical defatting and alkaline treatment process steps - Pseudorabies Virus (PrV). Eurofins Lancaster Labs, 2012.
- Study # NG-1292705: Virus clearance study for porcine peritoneal tissue for the chemical defatting and alkaline treatment process steps - Reovirus (Reo-3). Eurofins Lancaster Labs, 2012.
- Study # NG-1292706: Virus clearance study for porcine peritoneal tissue for the chemical defatting and alkaline treatment process steps - Porcine Parvovirus (PPV). Eurofins Lancaster Labs, 2012.
- Study # NG-1292709: Virus clearance study for porcine peritoneal tissue for the chemical defatting and alkaline treatment process steps - Bovine Viral Diarrhea Virus (BVDV). Eurofins Lancaster Labs, 2012.
- Study # NG-1282423: Virus clearance study for porcine peritoneal tissue for the gamma irradiation process - Bovine Viral Diarrhea Virus (BVDV). Eurofins Lancaster Labs, 2012.
- Study # NG-1282424: Virus clearance study for porcine peritoneal tissue for the gamma irradiation process - Porcine Parvovirus (PPV). Eurofins Lancaster Labs, 2012.
- Study # NG-1282425: Virus clearance study for porcine peritoneal tissue for the gamma irradiation process - Pseudorabies Virus (PrV). Eurofins Lancaster Labs, 2012.
- Study # NG-1282426: Virus clearance study for porcine peritoneal tissue for the gamma irradiation process - Reovirus (Reo-3). Eurofins Lancaster Labs, 2012.
- Study # NG-1282427: Virus clearance study for porcine peritoneal tissue for the gamma irradiation process - Xenotropic Murine Leukemia Virus (XMuLV). Eurofins Lancaster Labs, 2012.
- Study # K1/M12/16: Evaluation of the inactivation of murine leukemia virus, pseudorabies virus, bovine viral diarrhea virus, and porcine parvovirus by ethylene oxide treatment in the manufacturing process of Matricel collagen membrane. Charles River Biopharmaceutical Services GmbH, 2017.

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 3 of 20

- Study # K2/M12/16: Evaluation of the inactivation of porcine parvovirus by acid treatment in the manufacturing process of Matricel collagen membrane. Charles River Biopharmaceutical Services GmbH, 2017.
- Best SM, and Kerr PJ, 2000. Coevolution of host and virus: The pathogenesis of virulent and attenuated strains of myxoma virus in resistant and susceptible European rabbits. *Virology* 267: 36-48
- Brock KV, Grooms DL, Ridpath J, Bolin SR, 1998. Changes in levels of viremia in cattle persistently infected with bovine viral diarrhea virus. *J Vet Diagn Invest* 10: 22-26
- Charlton KM, Casey GA, Webster WA, 1984. Rabies virus in the salivary glands and nasal mucosa of naturally infected skunks. *Can J Comp Med* 48: 338-339
- Dogherty RM, Di Stefano H, 1971. Sites of avian leukosis virus multiplication in congenitally infected chickens. *Cancer Res* 27: 322-332
- Ferm VH, and Kilham L, 1963. Mumps virus infection of the pregnant hamster. *Embryol. Exp. Morph* 11: 659-665
- Gailunas P, and Cottrill GE, 1966. Presence and persistence of foot-and-mouth disease virus in bovine skin. *J. Bacteriol.* 91: 2333-2338
- Gibbs, S, Ellis AE, Mead DG, Allison AB, Moulton JK, Howerth EW, Stallknecht DE, 2005. West Nile virus detection in the organs of naturally infected Blue Jays (*Cyanocitta cristata*). *J Wildlife Dis* 41: 354-362
- Kiasari BA, Fallah FO, Monavari SH, 2013. Detection of viruses in the peritoneum of the Iranian patients on peritoneal dialysis. *Int J Med Med Sci* 5: 529-534
- Lackner C, Leydold SM, Modrow J, Farcet MR, Grillberger L, Schaefer B, Anderle H, Kreil TR. Reduction of spiked porcine circovirus during the manufacture of a Vero cell-derived vaccine, *Vaccine* 32: 2056-2061 (2014)
- Melnick JL, and Wallis C. Problems related to the use of serum and trypsin in the growth of monkey kidney cells. *Dev. biol. Stand.* 37: 77-82 (1977)
- Markus-Sekura C, Richardson JC, Harston RK, Sane N, Sheets RL: Evaluation of the human host range of bovine and porcine viruses that may contaminate bovine serum and porcine trypsin used in the manufacture of biological products. *Biologicals* 2011, 39:359-369
- Meng XJ: Emerging and re-emerging swine viruses. *Transboundary and emerging diseases* 59 (Suppl.1):85-102 (2012).
- Pedersen NC, 2014. An update on feline infectious peritonitis: Diagnosis and therapeutics. *Vet J* 201: 133-141
- Potts BJ and Kivens WJ: Detecting and removing porcine viruses to prevent contamination of biopharmaceutical raw materials. *BioPharm May* 2000:26-31.
- Royer RL, Nawagitgul P, Halbur PG et al.: Susceptibility of porcine circovirus type 2 to commercial and laboratory disinfectants. *J. Swine Health and Production* 9: 281-284 (2001)
- Tompkins SM, Lo CY, Tumpey TM, Epstein SL, 2004. Protection against lethal influenza virus challenge by RNA interference in vivo. *PNAS* 101: 8682-8686
- US-Pharmacopoeia - National Formulary (*USP-NF (1072) Disinfectants and Antiseptics*)
- Van Alstine WG, Kanitz CL, Stevenson GW, 1993. Time and temperature survivability of PRRS virus in serum and tissues. *J Vet Diagn Invest* 5:621-622
- Zwettler D, Fink M, Revilla-Fernández S, Steinrigl A, Winter P, Koefer J, 2012. First detection of Hepatitis E virus in Austrian pigs by RT-qPCR. *Berl Münch Tierärztl Wochenschr* 125: 281-289 (DOI 10.2376/0005-9366-125-281)

4 METHODOLOGY OF THE LITERATURE REVIEW

A literature review in accordance with Annex A of ISO 22442-3 was performed by Dr Gregersen in order to analyse and evaluate the data on the elimination and/or inactivation of viruses for the risk analysis of the viral safety of the Matricel Dental Barrier Membrane.

	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 4 of 20

4.1 Objective

The objective of the literature review was two-fold.

1. To compile a list with porcine viruses by virus families and virus types that might be expected as potential viral contamination of the peritoneal membranes collected. All cross-references need to be provided to allow a critical evaluation by regulatory authorities.
2. To assess the potential viral contamination levels in porcine raw materials in general and for peritoneal membranes in detail. The assessment should conclude with a realistic worst-case estimate for an unnoticed viral contamination level of peritoneal membranes in healthy pigs. All cross-references need to be provided to allow a critical evaluation by regulatory authorities.

4.2 Identification of data

Data was taken from scientific publications, e.g. valid scientific evidence, scientific literature and authoritative publications. A systematic search of the literature was conducted in order to reduce the risk of introducing bias.

The literature search was performed in the time period of June 18 till July 10, 2016. The following search strategy was used:

For the list of porcine viruses, a publication that is based upon a very broad literature search (Markus-Sekura et al. 2011) was taken as a base. This paper is commonly referred to as a reliable source for a list of potential porcine and bovine viruses. Due to certain restrictions of this paper's selection, further viruses of porcine origin were added to the list, which were taken from textbooks and specific publications (see details and references given in Section 6 below.) Likewise, taxonomic and morphological criteria were extracted from standard textbooks as referenced in Section 6 below.

To assess potential viral contamination levels in tissues, a systematic literature search was conducted using the "NLM Gateway" of the US national Library of Medicines, which includes, but is not limited to the Medline data bank (<http://gateway.nlm.nih.gov>). Neither publication date limit, nor a language or subset selection was applied. However, "Consumer Health" was excluded as a category. In an NLM Gateway search the search terms are automatically combined together (AND as default Boolean operator). Unqualified terms typed into the query box are automatically mapped to the according MeSH (medical subject headings) terms and searched as such, while the original term entered is also searched as text words. Therefore, simple terms with a broader scope and no Boolean operators were used as preferred search entries. Different search terms and combinations were applied (e.g. virus titer/tissue, virus titer/organ (or specific predilection organs names, such as liver, lung, tonsils, spleen, thymus, kidney, blood, serum), organ/titer, virus isolation/titer, virus isolation/pig (or other commonly studied species, such as bovine, mouse, chicken)) and tried to obtain suitable results. Hits were evaluated for their potential relevance by checking titles and abstracts. Where details seem to indicate a positive result but remained unclear, the full article was checked.

4.3 Assessment and critical evaluation

The assessment and critical evaluation of the data obtained during the literature search are provided in sections 5 and 6 below. The selection of the relevant/model viruses is described for the four experimental studies in sections 8.3, 9.3, 10.3 and 11.3.

5 POTENTIAL VIRAL CONTAMINANTS OF PORCINE ORIGIN

A wide range of viruses may be found in pigs. Table 1 lists viruses that may be found in pigs or which might grow in porcine cells and have a human host range. The list was assembled from the following sources:

- Markus-Sekura et al. (2011): Evaluation of the human host range of bovine and porcine viruses that may contaminate bovine serum and porcine trypsin used in the manufacture of biological products. Biologicals 2011, 39:359-369, and

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 5 of 20

- Potts BJ and Kivens WJ: Detecting and removing porcine viruses to prevent contamination of biopharmaceutical raw materials. BioPharm May 2000:26-31, and
- Meng XJ: Emerging and re-emerging swine viruses. Transboundary and emerging diseases 2012, 59 (Suppl.1):85-102.

In addition, the table was completed by basic viral characteristics which were extracted from Mahy BWJ (ed.): A dictionary of virology, 3rd. edition, Academic Press, 2001 and other virology textbooks.

Table 1: Potential viral contaminants of porcine origin

Virus Family	Virus	Genome	Envelope	Size (nm)	Physico-chemical resistance
Adenoviridae	Porcine adenovirus	dsDNA	no	65-80	medium
Anelloviridae	Torque teno virus (TTV)	ssDNA	no	17-20	very high
Arteriviridae	Porcine respiratory and reproductive syndrome virus (PRRS-irus)	ssRNA	yes	50-65	low
Asfarviridae	African swine fever virus	dsDNA	yes	200	low
Astroviridae	Porcine astrovirus	ssRNA	no	28-30	medium-high
Bunyaviridae	Crimean Kongo fever virus, Hantaan virus Jamestown canyon virus, La Crosse virus	ssRNA	yes	80-120	low
Caliciviridae	Porcine enteric calicivirus, norovirus, sapovirus, San Miguel sea lion virus	ssRNA	no	35-40	medium
Circoviridae	Porcine circovirus 1, 2 (PCV 1, 2)	ssRNA	no	17	very high
Coronaviridae	Transmissible gastroenteritis virus, porcine epidemic diarrhea virus, porcine haemagglutinating encephalitis virus, porcine respiratory coronavirus, porcine torovirus, SARS coronavirus, bovine coronavirus	ssRNA	yes	75-160	low
Filoviridae	Ebola reston virus	ssRNA	yes	80 x 800 -1000	low
Flaviviridae	Hog cholera virus, bovine viral diarrhea virus (BVDV), Bungowannah virus, dengue virus, Ilheus virus, Japanese encephalitis virus, louping ill virus, Murray Valley virus, Powassan virus, Wesselbron disease virus, West Nile virus	ssRNA	yes	50-70	low
Hepeviridae	Hepatitis E virus	ssRNA	no	32-34	medium-high
Herpesviridae	Pseudorabies virus (PSRV), porcine cytomegalovirus, porcine lymphotropic virus, bovine herpesvirus 1	dsDNA	yes	120-200	medium
Orthomyxoviridae	Influenzavirus	ssRNA	yes	80-120	low
Paramyxoviridae	Parainfluenza virus, Menangle virus, Nipah virus, peste-des-petit-ruminants virus, rinderpest virus, Tioman virus	ssRNA	yes	150-200 or more	low
Parvoviridae	Porcine parvovirus (PPV), bocavirus and hokovirus	ssDNA	no	18-24	very high

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 6 of 20

Virus Family	Virus	Genome	Envelope	Size (nm)	Physico-chemical resistance
Picornaviridae	Encephalomyocarditis virus, foot-and-mouth disease virus, porcine enterovirus 1 (Teschén/ Talfan disease) porcine enterovirus 2-9, swine vesicular disease virus, Seneca Valley virus, Kubovirus	ssRNA	no	25-30	medium-high
Poxviridae	Swine pox virus	dsDNA	yes	250 x 450	medium-high
Reoviridae	Reovirus, porcine rotavirus, banavirus	dsRNA	no	65-75	high
Retroviridae	Porcine endogenous retrovirus	ssRNA	yes	80-120	low
Rhabdoviridae	Rabies virus, vesicular stomatitis virus	ssRNA	yes	70 x 175	low
Togaviridae	Eastern/Venezuelan encephalitis virus, Getah virus, Ross River virus	ssRNA	yes	50-70	low

The publication by Markus-Sekura et al. (2011) considers not only viruses of porcine origin but also viruses that might have been picked up as contaminants of biomedicinal (cell culture) processes. Thus, Table 1 includes various viruses which are not found in pigs but may grow in porcine cell cultures and also have a human host range, for example various flaviviruses or Ebola virus. On the other hand, the list assembled by Markus-Sekura et al. (2011) did not mention certain porcine viruses that do not have a human host range. Using the other sources mentioned above, these missing porcine viruses were added to Table 1 and include porcine adenovirus, arterivirus, astroviruses, bocavirus and hog cholera virus.

Several viruses listed in Table 1 do not occur in Germany, such as African swine fever virus and most of the listed paramyxoviridae and togaviridae. Other porcine viruses listed above do not or do no longer occur in Germany as a result of disease control measures. This includes foot and mouth disease virus (no disease notification since 1990), stomatitis vesicularis, swine vesicular disease (no reported case since 1990), hog cholera (latest reported cases in domestic pigs until 2006, in wild boar until 2010). Pseudorabies virus still exists in Germany with very low prevalence (15 cases reported 2000 - 2013). [Disease notification data extracted from the German „TierSeuchenInformationssystem „TSIS“; latest reported status from Sept.17, 2014 under <http://tsis.fli.bund.de/Reports/YearOverview.aspx>]. According to the OIE (Office International des Epizooties) criteria, Germany is also free of classical mammalian rabies virus infections since no such cases have been observed since February 2006. However, isolated cases of bat lyssavirus in bats have been recorded in more recent years. Furthermore, importation of mammalian rabies, for example via imported dogs, cannot be fully excluded.

Viruses which might still be found in domestic pigs in Germany include porcine members of the adenoviridae (porcine adenovirus), astroviridae (porcine astrovirus) anelloviridae, (torque teno sus virus, TTV), caliciviridae (porcine enteric calicivirus, noro- and sapovirus), circoviridae (porcine circovirus 1 and 2), coronaviridae (several porcine viruses), hepeviridae (hepatitis E virus, HEV), orthomyxoviridae (influenza virus), paramyxoviridae (parainfluenzavirus), parvoviridae (porcine parvo-, boca-, and hokovirus), picornaviridae (porcine enteroviruses), reoviridae (reo-, rota, and banavirus), and retroviridae (porcine endogenous retrovirus).

Several of these viral agents cause only mild transient symptoms or silent or chronic infections without visible disease symptoms and might thus be overlooked. This applies particularly to those viruses that have only recently been identified or are not well studied, such as astroviruses, TTV, circoviruses, HEV, parvoviruses, reovirus, rotavirus, banavirus and retroviruses. Some of these might become human pathogens (for example HEV, influenza virus, parainfluenza virus, banavirus), others are considered strictly species-specific or apathogenic for humans (for example porcine astroviruses, TTV, PCV1/2, porcine rotavirus). But irrespective of their human pathogenicity, for a medicinal product any of these viruses are unwanted contaminants, even if they are only considered as „impurities“.

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 7 of 20

For a viral safety risk assessment of a product of porcine origin, the zoonotic viruses and human pathogens, in particular HEV virus, need specific attention. Furthermore, stable viruses, such as reo-, parvo- and circoviruses should be considered in viral safety risk assessments, in particular since these viruses have also been experienced as contaminants of biomedicinal processes. Parvo- and circoviruses are very small, non-enveloped viruses with a single-stranded genome. They are extremely stable under various physicochemical conditions. Parvoviruses are commonly considered worst-case model viruses for virus inactivation/reduction studies. Reoviruses are non-enveloped, double-stranded RNA viruses with a double-layered capsid. They show a rather high stability against certain virus inactivation agents that specifically react with nucleic acids. Thus reoviruses are also good and commonly applied models for viral safety assessments. The same applies to pseudorabies virus, a double-stranded DNA virus, which may also show a fairly high resistance, in particular to nucleic acid destroying agents.

Fairly stable enveloped virus, such as pseudorabies virus or BVDV can be considered suitable model viruses for the wide range of enveloped viruses, as most other enveloped viruses are of lower physicochemical resistance (see Table 1). Reoviruses are good models for non-enveloped virus of medium to high stability. Parvoviruses are among the smallest known viruses and with their extremely high resistance they can be considered model viruses, which, when successfully inactivated or reduced, might even justify claims for inactivation or reduction of any non-enveloped virus or even the entire spectrum of known viruses. However, it is not advisable to use only a parvovirus for model studies because only few processes have the capacity to inactivate or reduce parvovirus to the desired degree.

6 POTENTIAL VIRAL CONTAMINATION LEVELS IN PORCINE RAW MATERIALS

For acute virus infections with a specific and high tissue tropism and in organs of highest virus concentrations maximum virus titers were reported to amount up to 9 log₁₀ infectious units per gram tissue or per mL tissue extract (see for example Best and Kerr (2000) for myxoma virus in rabbits, or Gibbs et al. (2005) for West Nile virus in birds). Such high titers seem to mark a maximum that is rarely found and may occur only in specific organ samples and during a short period of time during the most acute phase of disease. Even for most viral predilection sites, average values were much lower.

Typical high-end titers in predilection tissues during the acute phase of a cyclic viral infection are normally in the range of 5 to 7 log₁₀ infectious units /mL or gram tissue extract (see for example Ferm and Kilham (1963) for mumps virus in hamsters, Dougherty and Di Stefano (1967) for avian leukosis virus in chicken, Charlton et al. (1984) for rabies virus in skunks, Tompkins et al. (2004) for mice challenged with influenza virus, Gibbs et al. 2005 for West Nile virus in various organ tissues).

Likewise, for viruses isolated from infected bovine or porcine specimens, virus titers were also in the range of up to 7 log₁₀ infectious units /mL or gram tissue. Brock et al. (1998) studied the levels of bovine viral diarrhea virus in serum collected from 4 persistently infected calves at birth until 70 days of age and found virus titers between 1.0 and 7.0 cell culture infective doses₅₀/ml. Gailunas and Cottrill (1966) studied the presence and persistence of foot-and-mouth disease virus in bovine skin and reported highest average virus titers of 3.6 log₁₀ plaque-forming units (PFU) per g of skin. Highest virus titers were found in blood samples in very early phases of the infection with titers up to 5.5 log₁₀ PFU/ml. Van Alstine et al. (1993) isolated PRRS virus from infected pigs and reported average virus titers of 5 log₁₀/ml for lung tissue (range 3 - 6), 2.8 log₁₀ for spleen and thymus tissues (range <2 - 4), and 2.6 log₁₀ for serum samples (range 2 - 3).

Using a quantitative PCR, Zwettler et al. (2012) studied the hepatitis E virus load of different virus strains in organs and samples of infected pigs. The PCR method applied did not differentiate between free viral nucleic acids and particle associated genomes and thus might overrate the true viral particle load and infectivity. The authors reported maximum HEV copy numbers of 6.2 log₁₀ for bile and liver and 4.1 log₁₀ for kidney samples. The virus loads of the positive reference samples, a bile and a fecal sample, were 5.7 and 5 log₁₀, respectively.

Compared to the viral loads quoted above, the potential viral load of peritoneal membranes collected from healthy pigs after veterinary in vivo and post mortem inspection is most likely at the very low end. The peritoneal membrane (lining of the abdominal cavity) is collected from healthy pigs and infections of peritoneal membranes would be clearly visible. Furthermore, peritoneal membranes are not a predilection site for viral infections. In

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 8 of 20

pathological, histological or immunohistochemical studies, signs of infection or viral antigens can normally be detected in the parenchyma of the affected abdominal organs or around blood vessels in interstitial tissues, whereas the organ-surrounding peritoneal membranes are normally not affected. There is only one virus that regularly causes peritonitis: the feline infectious peritonitis coronavirus. However, the peritonitis symptoms of this disease seem to be mainly caused by the deposition of inflammatory exudates and fibrinous immune complexes on peritoneal membranes, while the virus is found in various tissues and tends to replicate preferably in macrophages (Pedersen 2014). While cases of peritonitis are mostly associated with bacterial infections, peritoneal membranes may nevertheless also be "contaminated" by viruses during a systemic viral infection. For example, cytomegalovirus and BK-polyomavirus, viruses known to cause chronic disease with distribution of the virus in various organs and body fluids, have been found in cases of peritonitis in peritoneal dialysis patients (Kiasari et al. 2013).

Based upon the accumulated data for viral loads of acute or chronically infected tissues, a viral contamination level $7 \log_{10}$ infectious units per gram tissue appears to be a worst-case assumption for susceptible porcine tissues infected by an acute virus infection. Considering the specifics for peritoneal membranes of healthy pigs, this value might be too high and $6 \log_{10}$ might be a more appropriate, realistic worst-case estimate for an unnoticed viral contamination level of peritoneal membranes of healthy pigs.

7 VIRAL SAFETY MEASURES

Apart from manufacturing under controlled conditions in a clean room environment by healthy employees to control the contamination of the product in the Matricel facility, the viral safety of the DBM products is based on two pillars:

The first line of defence in the control of viral contamination of the products is the harvesting of porcine tissues under strict veterinary control in compliance with ISO 22442-2 to prevent tissue harvesting from diseased animals that are potentially infected with viruses.

The second line of defence is the incorporation of viral inactivating sub-processes within the product manufacturing process. The Matricel manufacturing/purification process of porcine peritoneal tissues inherits different process elements with the potential for viral inactivation.

These are:

- Sub-process 1: For disinfection a treatment of the tissues with at least 1.0 N NaOH and 4.0 M NaCl is performed at least at 20°C for at least 90 minutes.
- Sub-process 2: For chemical defatting the tissues are treated in 98% to 100% acetone for at least 18 hours at ambient temperature.
- Sub-process 3: During the alkali treatment the tissues are purified with at least 0.1 N NaOH (pH > 13) for 22 to 25 hours at ambient temperature, and
- Sub-process 4: Subsequently an acidic treatment with at least 0.1 N HCl (pH < 1.0) is performed for 20 to 22 hours at ambient temperature.
- Sub-process 5: Terminal sterilization by ethylene oxide (EO) gas treatment (HA2 routine cycle 05; Steam injection 16-25mbar, EO concentration 644-686 mg/l, treatment time 239-241 min., Treatment temperature 42-50°C).

Viral reduction/inactivation studies have been performed for sub-processes 2, 3, 4, and 5. The experimental data will be summarized in section 8, 9, 10 and 11 below. The results will be discussed and a critical evaluation of the overall residual risk acceptability will be provided in section 12 of this Risk Analysis.

	Risk Analysis	Document: QP-RA-03005-B
		Release: 28.Feb.17
	Viral safety of Dental Barrier Membrane	Page: 9 of 20

8 VIRAL SAFETY DATA AVAILABLE FOR BM COLLAGEN MANUFACTURING/PURIFICATION PERFORMED BY NEULAB (NL-STUDY)

8.1 Background information

The inactivation of the viruses by the acetone treatment (sub-process 2 described above) and NaOH treatment (sub-process 3 described above) was studied in co-operation with NewLab BioQuality AG (Gottfried-Hagen-Str. 62, D-51105 Köln, Germany). All sub-process steps were adjusted to worst case conditions (e.g. shortest treatment time, lowest concentration of chemicals and temperature). Test items were spiked with model viruses and before and after the respective treatments, samples were withdrawn and assayed for their virus titer by endpoint titration.

The study was conducted according to the following guidance documents valid in 2005:

- Anforderungen an Validierungsstudien zum Nachweis der Virussicherheit von Arzneimitteln aus menschlichem Blut oder Plasma (Paul Ehrlich Institut und Bundesgesundheitsamt; May 4, 1994)
- Note for Guidance on Virus Validation Studies (CPMP/BWP/268/95, Final Version 2; February 29, 1996)
- Note for Guidance on Plasma Derived Medicinal Products (CPMP/BWP/269/95, Revision 3, January 25, 2001)
- EN 12442-3:2000 (Animal tissues and their derivatives utilized in the manufacture of medical devices – Part 3: Validation of the elimination and/or inactivation of viruses and transmissible agents).

8.2 Selection of viruses

In discussions with the viral experts of NewLab and according to the Note for Guidance on Virus Validation Studies the following selection criteria were considered during selection of the viruses:

1. Host of the virus (species);
2. Genome of the virus (DNA or RNA) and;
3. Structure of the virus (enveloped or non-enveloped).

The selected viruses should be preferably of porcine origin, i.e. the same origin as the source tissue from which the collagen is isolated. Secondly, the viral genome should include a double stranded virus. Thirdly, due to significant differences in resistance between enveloped and non-enveloped viruses, a non-enveloped and a non-enveloped virus should be studied. Based upon those criteria, pseudorabies virus and porcine parvovirus were chosen as model viruses.

Pseudorabies virus (PRV, suid herpesvirus 1). This virus belongs to the herpesviridae family, alphaherpesvirinae subfamily, varicellovirus genus. The virions are enveloped, have icosahedral symmetry with 120 to 200 nm diameter. The genome consists of a double-stranded DNA. This virus has a medium resistance to various physicochemical treatments.

Porcine parvovirus (PPV). This virus belongs to the parvoviridae family, parvovirinae subfamily, parvovirus genus. The virions are icosahedral non-enveloped particles with a 18-26 nm diameter. The genome consists of a single stranded DNA. This virus has a very high resistance to various physicochemical treatments.

8.3 Definition of worst case process and experimental conditions

Chemical defatting - acetone treatment: BM collagen manufacturing involves the treatment of porcine tissue with acetone, as part of the chemical defatting step in the purification process. The chemical defatting step consists of three successive acetone washes (Step #1, Step #2, and Step #3) and a water wash. The process step that was validated is Step #2 of the defatting sub-process. Step #2 of this sub-process describes that the porcine tissue is treated at room temperature in 15 liters of acetone for a minimum of 18 hours and a maximum of 14 days under continuous stirring during working hours. The density of the acetone should be below 0.805 g/cm³.

matricel	Risk Analysis	Document: QP-RA-03005-B
		Release: 28.Feb.17
	Viral safety of Dental Barrier Membrane	Page: 10 of 20

Based on the description above, worst case conditions for the acetone treatment were defined as follows: Treatment with acetone (density not below 0.805 g/cm³) for 18 hours at 18 °C.

NaOH treatment: BM collagen manufacturing involves the treatment of porcine tissue with NaOH. The sub-process step to be validated describes that the porcine tissue is treated at room temperature in 25 liters of a NaOH solution for a minimum of 22 hours and a maximum of 25 hours. The NaOH solution should have a concentration equal to or more than 0.1 N.

Based on the description above, worst case conditions for the 0.1 N NaOH treatment are defined as follows: Treatment with a NaOH solution (concentration 0.1 N) for 22 hours at 18 °C.

Justification for the scaled down parameters and details on the loading of the test item with virus particles can be found in the NewLab report. The critical parameters defined above are monitored during routine production. This guarantees that the experimental conditions during viral inactivation testing are also met during routine manufacturing/purification conditions.

8.4 Calculations

In the NL-study, test items were spiked with model viruses at a defined titer. Before and after the respective treatments, samples were withdrawn and monitored for virus titer by endpoint titration. The reduction factors achieved by the respective treatments were calculated as follows:

$$\text{Reduction factor } R = \log_{10} A_0 - \log_{10} A_n$$

A_0 : total viral load in the control (media/NaCl)

A_n : total viral load in the process sample

The results reported here are the worst-case experimentally estimated reduction factors obtained including the 95% confidence limits. A detailed description of all calculations can be found in the NewLab report.

8.5 Results

All viral clearance studies were performed in duplicate. A summary of the results is presented in Table 2 below. A detailed description of the results of the study can be found in the NewLab report.

Chemical defatting - acetone treatment: PRV was inactivated to the limit of detection by acetone treatment. PPV was not inactivated by acetone treatment. The viral input could almost completely be recovered from the spiked test item.

NaOH treatment: Both, PRV as well as PPV, were inactivated by 0,1 N NaOH treatment.

Table 2: Worst-case reduction factors ($\pm$ the 95% confidence limits) obtained in the NL-study.

Treatment	PRV		PPV	
	Run 1	Run 2	Run 1	Run 2
Acetone treatment	$\geq 4.61 \pm 0.28$	$\geq 4.55 \pm 0.20$	0	0
NaOH treatment	3.97 ± 0.29	3.97 ± 0.22	6.00 ± 0.26	4.69 ± 0.29

	Risk Analysis	Document: QP-RA-03005-B
		Release: 28.Feb.17
	Viral safety of Dental Barrier Membrane	Page: 11 of 20

9 VIRAL SAFETY DATA AVAILABLE FOR BM COLLAGEN MANUFACTURING/PURIFICATION PERFORMED BY LANCASTER LABS (LL-STUDY)

9.1 Background information

A second experimental viral inactivation study was performed in co-operation with Genzyme - a Sanofi Company, 64 Sidney Street, Cambridge, MA 02139, as a sponsor. Genzyme is a long-term partner of Matricel, and both companies have identical interest in establishing the viral inactivation capability of the manufacturing/purification process used by Matricel to purify porcine peritoneal tissues. The facility conducting the study was Eurofins Lancaster Laboratories, Incorporated, Pharmaceutical Science Division, Viral Safety and Clearance Services, 2425 New Holland Pike, Lancaster, PA 17601.

The LL-study expanded upon the NL-study described above, with the identical sub-processes selected for analysis: Chemical defatting using acetone treatment and alkaline treatment. Other than the NL-study, the LL-study investigated multiple wash steps within the chemical defatting process and multiple time points during the alkaline treatment. In addition, the LL-study further expands upon the Matricel study by testing inactivation of a panel of five viruses that vary in size, envelope, genome type, and expected resistance to chemical inactivation.

The LL-study was conducted in accordance with:

- Q5A viral safety evaluation of biotechnology products derived from cell lines of human or animal origin. Federal Register 1998 Sept 24;63(185):51075-51081.
- Note for Guidance on Virus Validation Studies: The design, contribution and interpretation of studies validating the inactivation and removal of viruses. European Agency for the Evaluation of Medicinal Products, CPMP/BWP/268/95, February 14, 1996.
- Quality of biotechnological products: viral safety evaluation of biotechnological products derived from cell lines of human or animal origin (ICH Topic Q5A, Step 3 document) (CPMP/ICH/295/95).
- Medical devices utilizing animal tissues and their derivatives—Part 3: Validation of the elimination and/or inactivation of viruses and transmissible spongiform encephalopathy (TSE) agents (ANSI/AAMI/ISO 22442-3:2007).

9.2 Selection of viruses

The porcine origin of the membrane guided the selection of the model viruses for the clearance studies (Table 3). Pseudorabies virus (PRV) and porcine parvovirus (PPV) are of porcine origin. Porcine endogenous retrovirus (PERV) is a known gamma retrovirus found in swine; however, there is no validated infectivity assay for this virus. Xenotropic murine leukemia virus was selected as a representative gamma retrovirus for use as a model for PERV. Reovirus and bovine viral diarrhea virus were selected as non-specific model RNA viruses. These additional viruses were tested as models for possible contaminants arising during production. All assays used in this study have been validated according to ICH Q2.

Table 3: Viruses selected for the LL- study.

Virus	Family	Genome	Envelope	Size (nm)	Physicochemical resistance
Xenotropic murine leukemia virus (XMuLV)	Retro	RNA	Yes	80-110	Low
Pseudorabies virus (PRV)	Herpes	DNA	Yes	120-200	Medium
Reovirus type 3 (Reo-3)	Reo	RNA	No	60-80	High
Porcine parvovirus (PPV)	Parvo	DNA	No	18-24	Very high
Bovine Viral Diarrhea Virus (BVDV)	Flavi	RNA	Yes	50-70	Low

	Risk Analysis	Document: QP-RA-03005-B
		Release: 28.Feb.17
	Viral safety of Dental Barrier Membrane	Page: 12 of 20

9.3 Definition of worst case process and experimental conditions

Chemical defatting - acetone treatment: The chemical defatting sub-process consists of three successive acetone washes (Step #1, Step #2, and Step #3) and a water wash. The moisture content of the acetone was expected to be the key parameter in the inactivation process. Historical values of acetone wash density (related to moisture content) were used to determine the worst-case conditions for each wash. Three standard deviations above the average moisture content were used as a basis for Steps #1 and #3. Step #2 is required by protocol not to exceed 0.805 g/cm³. The shortest allowable process duration, the lowest allowable wash temperature and the mildest mixing regime possible were used as other worst case approximations. A summary of the experimental conditions is provided below.

Step #1: After spiking and incubation, all six test samples were placed in conical tubes and washed with 17.4 % v/v moisture content in acetone for 90 minutes. Two samples were removed for viral load testing.

Step #2. The four samples that were not removed for analysis were transferred to fresh tubes and gently washed in 5% v/v moisture content in acetone for 2-2.5 hours, then transferred to an 18°C water bath for the balance of 18 hours. The four tubes were then transferred to a 33°C water bath for 30 minutes. Two samples were removed for viral load testing.

Step #3. The two remaining samples were transferred to fresh tubes and washed in 1% v/v moisture content in acetone for 90 minutes. The two tubes were transferred to a 33°C water bath for 30 minutes. The two samples were then removed for viral load testing.

NaOH treatment: The alkaline treatment sub-process consists of a 0.1 NaOH wash followed by HCl neutralization and sodium bicarbonate buffering. The shortest allowable process duration, the lowest allowable wash temperature, the smallest allowable solution volume and the mildest mixing regime possible were used as other worst case approximations. A summary of the experimental conditions is provided below.

After spiking and incubation, the six test samples were placed in conical tubes and 0.1N NaOH was added. The tubes were placed on a tube rotator for 1h, 5h, and 22h. At each time point, two samples were removed and immediately buffered with 0.357M sodium bicarbonate, neutralized with 1N HCl, and submitted for viral activity analysis.

Justification for the scaled down parameters and details on the loading of the test item with virus particles can be found in the LL-study reports. The critical parameters defined above are monitored during routine production. This guarantees that the experimental conditions during viral inactivation testing are also met during routine manufacturing/purification conditions.

9.4 Calculations

Virus titers and viral clearance factors were calculated according Eurofins Lancaster Laboratories procedure. The virus clearance factor (VCF) for each step was calculated as follows:

$$VCF = \log_{10} \left(\frac{(\text{Volume}) \times (\text{titer before processing})}{(\text{Volume}) \times (\text{titer after processing})} \right)$$

The Start sample titer was used for the "before processing" value as recommended in Q5A Guidance and is reported here in this risk analysis.

Additional VCFs for supernatants "after processing" and VCFs using hold samples or shipping controls as the comparator were calculated. These values are reported in the applicable original reports obtained from Lancaster Labs.

9.5 Results

All viral clearance studies were performed in duplicate. A summary of the results is presented in Table 4 and Table 5 below. A detailed description of the results can be found in the LL-study reports.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 13 of 20

Table 4: Viral clearance factors of the chemical defatting sub-process, using start samples as the comparator.

Test sample	Run 1		Run 2	
	VCF	95% Confidence Limit	VCF	95% Confidence Limit
XMuLV Step #1	≥5.22	0.07	≥5.07	0.13
XMuLV Step #2	≥5.22	0.07	≥5.07	0.13
XMuLV Step #3	≥5.18	0.07	≥4.98	0.13
PRV Step #1	≥5.71	0.03	≥5.93	0.06
PRV Step #2	≥5.71	0.03	≥5.93	0.06
PRV Step #3	≥5.71	0.03	≥5.93	0.06
Reo-3 Step #1	6.74	0.88	≥6.25	0.05
Reo-3 Step #2	≥6.26	0.09	≥6.25	0.05
Reo-3 Step #3	≥6.26	0.09	≥6.25	0.05
PPV Step #1	0.12	0.14	0.09	0.15
PPV Step #2	0.13	0.08	0.09	0.13
PPV Step #3	0.37	0.20	0.30	0.25
BVDV Step #1	≥5.09	0.11	≥5.13	0.08
BVDV Step #2	≥5.09	0.11	≥5.13	0.08
BVDV Step #3	≥5.09	0.11	≥5.13	0.08

Table 5: Tissue viral clearance factors of the alkaline treatment sub-process, using start samples as the comparator.

Test sample	Run 1		Run 2	
	VCF	95% Confidence Limit	VCF	95% Confidence Limit
XMuLV 1h	≥4.89	0.15	≥5.07	0.14
XMuLV 5h	≥4.89	0.15	≥5.07	0.14
XMuLV 22h	≥4.89	0.15	≥5.07	0.14
PRV 1h	≥5.89	0.09	≥5.71	0.12
PRV 5h	≥5.89	0.09	≥5.71	0.12
PRV 22h	≥5.89	0.09	≥5.71	0.12
Reo-3 1h	6.77	0.89	6.10	0.46
Reo-3 5h	≥6.29	0.14	6.80	0.88

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 14 of 20

Test sample	Run 1		Run 2	
	VCF	95% Confidence Limit	VCF	95% Confidence Limit
Reo-3 22h	6.47	0.62	≥6.32	0.13
PPV 1h	4.96	0.09	5.09	0.07
PPV 5h	4.81	0.14	4.94	0.06
PPV 22h	5.89	0.10	5.80	0.06
BVDV 1h	≥5.46	0.05	≥5.39	0.10
BVDV 5h	≥5.46	0.05	≥5.39	0.10
BVDV 22h	≥5.46	0.05	≥5.39	0.10

10 VIRAL SAFETY DATA FOR THE ACIDIC TREATMENT USED DURING BM COLLAGEN MANUFACTURING PERFORMED BY CHARLES RIVER (CRL-2 STUDY)

10.1 Background information

An additional experimental viral inactivation study was performed by Charles River Biopharmaceutical Services GmbH, Gottfried-Hagen-Str. 20, 51105 Köln, Germany (CRL2-study, study no. K2/M12/16). The CRL2-study investigated the viral clearance achieved by the acidic treatment.

The CRL2-study was conducted in accordance with:

- "Points to consider in the characterization of cell lines used to produce biologicals"
- "Note for guidance on virus validation studies"
- "Note for guidance on quality of biotechnological products" (ICH Q5A)
- "Points to consider in the manufacturing and testing of monoclonal antibody products for human use"
- "Guideline on virus safety evaluation of biotechnological investigational medicinal products"
- "Note for guidance on plasma derived medicinal products"
- Validation of the elimination and/or inactivation of viruses and transmissible spongiform encephalopathy (TSE) agents.

10.2 Selection of viruses

Porcine Parvovirus was chosen as a model virus for this study, as it represents extremely stable viruses with a high resistance against acid pH values. Furthermore, parvoviruses are considered potential porcine viral contaminants of the starting material.

Table 6: Viruses selected for the CRL2- study.

Virus	Family	Genome	Envelope	Size (nm)	Physicochemical resistance
Porcine parvovirus (PPV)	Parvo	DNA	No	18-24	Very high

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 15 of 20

10.3 Definition of worst case process and experimental conditions

The acidic treatment consisted of a 0.1N HCl treatment followed by direct recovery of the remaining viral load by a threefold wash with medium and an additional wash/virus elution step with 1 M NaCl. In control studies these conditions showed almost complete virus recovery. Virus reduction was determined after evaluating the worst case conditions of the manufacturing process. The shortest allowable process duration and the mildest mixing regime possible were used as other worst case approximations.

After spiking the test item and air drying, 0.1N HCl was added. The materials were mildly agitated using a stirrer for 1h, 5h, and 20h at room temperature. At each time point, two samples were removed, passed through the recovery procedure and tested for residual virus.

Justification for the scaled down parameters are on file at Matricel and details on the loading of the test item with virus particles can be found in the CRL2-study report. The critical parameters are monitored during routine production. This guarantees that the experimental conditions during viral inactivation testing are also met during routine manufacturing/purification conditions.

10.4 Calculations

For the calculation of the reduction factor of a process fraction the virus load of the process fraction is referred to the virus load of the load fraction (reduction factor).

$$R = \log_{10} A_0 - \log_{10} A_n$$

R: reduction factor
 A_0 : reference virus load
 A_n : sample / process virus load

Some samples were analyzed by large volume plating in addition to the endpoint titration. Because the tested volume in the LVP is much higher than the tested volume in the titration, the result of the LVP is statistically more reliable. Consequently, the reduction factor reported for such a sample is calculated using the viral load determined by LVP.

The load sample titer was used for the "reference virus load" value as recommended in Q5A Guidance and is reported here in this risk analysis.

10.5 Results

All viral clearance studies were performed in duplicate. A summary of the results is presented in Table 7 below. A detailed description of the results can be found in the CRL2-study report.

Table 7: Tissue viral clearance factors of the acidic treatment sub-process, using start samples as the comparator.

Test sample	Run 1		Run 2	
	LRF	95% Confidence Limit	LRF	95% Confidence Limit
PPV 1h	3.92	0.32	≥4.11	0.27
PPV 5h	≥6.55	0.30	6.77	0.27
PPV 20h	≥6.56	0.30	≥6.33	0.27

	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 16 of 20

11 VIRAL SAFETY DATA FOR THE EO STERILIZATION METHOD USED DURING DBM MANUFACTURING PERFORMED BY CHARLES RIVER (CRL-1 STUDY)

11.1 Background information

An additional experimental viral inactivation study was performed by Charles River Biopharmaceutical Services GmbH, Gottfried-Hagen-Str. 20, 51105 Köln, Germany (CRL1-study, study no. K1/M12/16) by studying the viral clearance achieved by the EO sterilization method using a similar set of viruses as in the LL-Study above (the Reo-3 virus was not included since for this virus already sufficient reduction was demonstrated).

The CRL1-study was conducted in accordance with:

- "Points to consider in the characterization of cell lines used to produce biologicals"
- "Note for guidance on virus validation studies"
- "Note for guidance on quality of biotechnological products" (ICH Q5A)
- "Points to consider in the manufacturing and testing of monoclonal antibody products for human use"
- "Guideline on virus safety evaluation of biotechnological investigational medicinal products"
- "Note for guidance on plasma derived medicinal products"
- Validation of the elimination and/or inactivation of viruses and transmissible spongiform encephalopathy (TSE) agents.

11.2 Selection of viruses

The porcine origin of the membrane guided the selection of the model viruses for the clearance studies, as already presented in section 9.2 above. Pseudorabies virus (PRV) and porcine parvovirus (PPV) are of porcine origin. Porcine endogenous retrovirus (PERV) is a known gamma retrovirus found in swine; however, there is no validated infectivity assay for this virus. Xenotropic murine leukemia virus was selected as a representative gamma retrovirus for use as a model for PERV. Bovine viral diarrhea virus was selected as a non-specific model RNA virus. These additional viruses were tested as non-specific models for possible contaminants arising during production. All assays used in this study have been validated according to ICH Q2.

Table 8: Viruses selected for the CRL1- study.

Virus	Family	Genome	Envelope	Size (nm)	Physicochemical resistance
Xenotropic murine leukemia virus (MuLV)	Retro	RNA	Yes	80-110	Low
Pseudorabies virus (PRV)	Herpes	DNA	Yes	120-200	Medium
Porcine parvovirus (PPV)	Parvo	DNA	No	18-24	Very high
Bovine Viral Diarrhea Virus (BVDV)	Flavi	RNA	Yes	50-70	Low

11.3 Definition of worst case process and experimental conditions

Routine EO sterilization is performed by HA2 Medizintechnik GmbH, Am Bahndamm 11, 38820 Halberstadt, Germany, using HA2 sterilization cycle 05, under mixed load conditions. The detailed description of cycle 05 is described in HA2 document FB-FG-109, current revision 07, released 14.Aug.2015. For process validation under worst case conditions, HA2 has defined half-cycle 05 conditions in accordance with the requirements of ISO 11135 (Annex B Conservative determination of lethal rate of the sterilization process – Overkill approach). These half-cycle 05 conditions have also been used as worst case conditions during the CRL1-study. A detailed description of the routine cycle 05 and half-cycle 05 conditions is provided in Matricel document QP-VR-03009-A. A summary is provided in Table 9 below.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 17 of 20

Table 9: Worst case conditions used during EO sterilization

Aspect of sterilization	Worst case conditions	Rationale
Mixed load configuration	1 pallet space of Sterilization Chamber 6 was used for the box containing the test items. The remaining 5 pallet spaces were filled with the worst case dummy pallets used by HA2 during validation activities.	The worst case dummy pallets contain 250 kg of PVC, a plastic known for its high EO absorbance. This created a worst case EO environment during the EO sterilization step.
Pre-stock holding	Was not performed.	This created worst case conditions when compared to routine cycle 05 conditions: T = $\geq 15^{\circ}\text{C}$ t = $\geq 12\text{h}$
Pre-conditioning	The "half-cycle 05 conditions" were used: T = 30.0-39.9°C RH = 0-50% t = 6-7h	This created worst case conditions when compared to routine cycle 05 conditions: T = 40.0-50.0°C RH = 50-85% t = 7-48h Note: This step is crucial to create the moist/warm conditions to allow the EO to become an active sterilizing agent. Therefore it cannot be deleted.
EO sterilization	The "half-cycle 05 conditions" were used: Steam = 0-16mbar EO = 626-651 mg/l t = 119-121 min. T = 32.0-41.9°C	This created worst case conditions when compared to routine cycle 05 conditions: Steam = 16-25mbar EO = 644-686 mg/l t = 239-241 min. T = 42-50°C
Post-conditioning	The "half-cycle 05 conditions" were used: T = 30.0-36.9°C t = 19-20h	This created worst case conditions when compared to routine cycle 05 conditions: T = 37.0-47.0°C t = 20-48h Note: This step ensured that residual EO is effectively removed from the test items, to prevent prolonged activity of EO during transportation and interference with the viral clearance study (toxicity for indicator/propagation cells)
Quarantine storage	Was not performed.	Is not performed during routine sterilization of Matricel products

Details on the loading of the test item with virus particles can be found in the CRL2-study report. The critical parameters defined above are monitored during routine production. This guarantees that the experimental conditions during viral inactivation testing are also met during routine manufacturing/purification conditions.

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

matricel	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 18 of 20

11.4 Calculations

For the calculation of the reduction factor of a process fraction the virus load of the process fraction is referred to the virus load of the load fraction (reduction factor).

$$R = \log_{10} A_0 - \log_{10} A_n$$

R: reduction factor
A₀: reference virus load
A_n: sample / process virus load

Some samples were analyzed by large volume plating in addition to the endpoint titration. Because the tested volume in the LVP is much higher than the tested volume in the titration, the result of the LVP is statistically more reliable. Consequently, the reduction factor reported for such a sample is calculated using the viral load determined by LVP.

The load sample titer was used for the "reference virus load" value as recommended in Q5A Guidance and is reported here in this risk analysis.

11.5 Results

All viral clearance studies were performed in duplicate. A summary of the results is presented in Table 10 below. A detailed description of the results can be found in the CRL1-study report.

Table 10: Tissue viral clearance factors of the EO sterilization process, using start samples as the comparator.

Test sample	Run 1		Run 2	
	LRF	95% Confidence Limit	LRF	95% Confidence Limit
MuLV	≥3.99	0.22	≥4.11	0.25
PRV	≥5.00	0.31	≥4.58	0.32
PPV	6.33	0.24	6.98	0.23
BVDV	≥5.24	0.28	≥4.94	0.27

12 DISCUSSION AND CRITICAL EVALUATION OF THE OVERALL RESIDUAL RISK ACCEPTABILITY

A summary of the virus inactivation data obtained during our studies along with the cumulative overall log viral clearance factor is given in Table 11 below.

The LL-study experimentally confirmed the results found in the NL-study. Compared to the NL-study, the LL-study investigated multiple wash steps within the chemical defatting process and multiple time points during the alkaline treatment.

The chemical defatting process step is highly effective in inactivating virus, with viral clearance factors of ≥ 4.9 log₁₀ for all of the viruses tested except for PPV, a non-enveloped virus. No XMuLV, PRV, Reo-3 or BVDV virus was detected after chemical defatting Step #1, despite its short duration (90 minutes), high moisture content (17.4 % v/v) and mild temperature. The two subsequent chemical defatting steps are of longer duration, lower moisture content, and include treatment at elevated temperature, all of which should contribute to further reducing the activity of viruses susceptible to acetone exposure.

Alkaline treatment proved to be a robust method for viral inactivation, i.e. significant viral clearance is already observed during the early stages (1-5 hours) of the treatment of 22 hours. Alkaline treatment is effective against

\\Sbs\Document Control\QP documents\QP-RA\QP-RA-03005\QP-RA-03005-B.docx

	Risk Analysis	Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane	Release: 28.Feb.17
		Page: 19 of 20

all five viruses, including the non-enveloped viruses. No XMuLV, PRV, or BVDV virus was detected at the end of this process step. The minimum viral clearance factor determined at the end of the process ranged from 3.97 to ≥ 6.32 for the five viruses tested.

Acidic treatment proved to be a robust method for viral inactivation i.e. significant viral clearance is already observed during the early stages (1 or 5 hours) of the treatment of 20 hours. Acidic treatment was tested only for PPV, the virus with highest physico-chemical resistance effective. No PPV virus was detected at the end of this process step. The minimum viral clearance factor determined at the end of the process was ≥ 6.33 .

EO sterilization was effective against all viruses tested (MuLV, PRV, PPV and BVDV). Under worst case conditions, no virus was found after the treatment for MuLV, PRV and BVDV, while a log reduction reduction of 6.33 was found for PPV.

Based upon the data summarized above, the DBM manufacturing process involves four proven and highly effective virus inactivation/reducing steps for enveloped viruses (MuLV, PRV, and BVBV), including those enveloped viruses that are more resistant (e.g. PRV) than the majority of the lipid containing viruses. The cumulative viral clearance factor for enveloped viruses was in any case greater than $13.10 \log_{10}$.

Reoviruses, mainly representing medium and highly resistant, non-enveloped viruses with a more complex architecture, such as adenoviruses, poxviruses and other members of the reoviridae family (rota- and banavirus) are strongly affected by all three tested virus reducing/inactivation steps, resulting in a cumulative viral clearance factor of at least $12.57 \log_{10}$. Reo-3 was included in only two of the four experimental studies, therefore this cumulative viral clearance factor of at least $12.57 \log_{10}$ must be considered an underestimation of the real viral clearance.

PPV, one of the most resistant viruses known and representative for any naked virus of high or very high resistance (e.g. anelloviruses, astroviruses, circoviruses, hepeviruses, picornaviruses, and other members of the parvoviridae family, namely boca- and hokoviruses) is effectively inactivated by three of the tested methods with a cumulative viral clearance factor of at least $17.35 \log_{10}$.

Compared to a realistic worst-case viral contamination level of $6 \log_{10}$ for the peritoneal membrane raw material, the cumulative viral clearance factors indicate a very high safety margin of ≥ 6 logs for all viruses including extremely stable, small, non-enveloped viruses. A safety margin of 6 logs means that even if a high viral contamination of the raw material should occur, the DBM process would reduce the contaminating virus to less than one millionth of an infectious dose.

We conclude that the DBM manufacturing process reproducibly and reliably renders a safe product regarding its viral safety.

matricel	Risk Analysis		Document: QP-RA-03005-B
	Viral safety of Dental Barrier Membrane		Release: 28.Feb.17
			Page: 20 of 20


Table 11: Cumulative viral clearance factors obtained in all experimental viral inactivation studies performed.

Treatment	XMuLV		PRV			Reo-3		PPV			BVDV	
	LL-study	CRL1-study	NL-study	LL-study	CRL1-study	LL-study	CRL1-study	NL-study	LL-study	CRL1 and CRL2-study	LL-study	CRL1-study
Chemical defatting - Acetone treatment	≥4.98		≥4.55	≥5.71		≥6.25		0*	0*		≥5.09	
NaOH treatment	≥4.89		3.97	≥5.71		≥6.32		≥4.69	5.80		≥5.39	
HCl treatment										≥6.33		
EO sterilization		≥3.99			≥4.58					6.33		≥4.94
Cumulated log clearance factor	≥13.86		≥13.10 - ≥16.00			≥12.57		≥17.35 - ≥18.46			≥15.42	

*values less than 1.00 were not included in the sum calculation
Only the lowest log reduction values of repeated runs of a study were included in the table.

matricel	Verification Report	Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1	Release: 27.Jan.14 Page: 1 of 19

PREPARED BY:

 L. Olde Damink

DATE: 27 Jan'14.

1. INTRODUCTION

The Matricel Dental Barrier Membrane (DBM) is a porcine membrane intended to be used in dental applications as a barrier membrane for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR) to support the growth of new tissue and bone. The verification study summarized in this report was designed to answer questions raised during the review of 510(k) premarket notification K123697 to FDA. The study was performed by Toxikon, 15 Wiggins Ave. Bedford MA, USA under number 13-02011-G1. The study conformed to the current FDA 21 CFR Part 58 Good Laboratory Practice for Nonclinical Laboratory Studies. This report is prepared to provide a summary of the final report obtained from Toxikon.

The purpose of the study was to evaluate the DBM membrane for its degradation rate after implantation over dental extraction sockets, compared to the predicate membrane Geistlich Resorbable Bilayer Membrane Bio-Gide. The study assessed tissue integration and reaction to the membranes. In addition the amount of bone formation at the adjacent implanted bone defect was measured in order to verify that the tissue reaction to the membranes did not interfere with the adjacent bone formation taking place at the grafted site. For this purpose the dental extraction sockets were filled with bone graft material (Bio-Oss® spongiuous bone substitute). Grafted sites were evaluated for bone regeneration, bone resorption and remodeling, ingrowth of surrounding tissues, inflammation, and effect on surrounding tissues.

2. REFERENCES

- Final GLP Report 13-02011-G1: Dental Barrier Membrane Implantation in Beagle Dogs.
- Giannobile W, Nevins M. Osteology Guidelines for Oral and Maxillofacial Regeneration. Preclinical Models for Translational Research. Osteology Foundation. Quintessence Publishing 2011.

3. MATERIALS

For this animal study, the Geistlich Resorbable Bilayer Membrane Bio-Gide was selected as the predicate device. Both the Matricel DBM (also referred to as test article) and predicate device (also referred to as control article) are derived from porcine tissue and share similar technological and biochemical characteristics. The Matricel DBM and Bio-Gide membrane contain Type I with Type III collagen and elastin fibers. Furthermore, the available sizes and weight per unit surface area (surface weight) of the Matricel DBM are within the specifications of Bio-Gide. Both the proposed membrane and Bio-Gide are not chemically crosslinked and are resorbable over time. Both the proposed and the predicate devices are intended for use for GBR and GTR during specific dental procedures involving implantology and periodontology. The following materials were used during the study:

Test Article: Matricel Dental Barrier Membrane (Matricel GmbH), Lot/Batch Number REM1520-1301, Use-by date: 2016-03.

Control Article: Resorbable Bilayer Membrane Bio-Gide (Geistlich Pharma), Lot/Batch Number 00120797 Use-by date: 2015-08.

Bone Graft Material: Bio-Oss® spongiuous bone substitute, large granules, 1-2 mm (Geistlich Pharma), Lot/Batch Number 00120782 Use-by date: 2015-11 and Lot/Batch Number 00130073 Use-by date: 2015-09.

matricel	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 2 of 19

4. RATIONALE FOR THE DESIGN OF THE STUDY

As part of the experiments to demonstrate the substantial equivalence in terms of safety and performance of the Matricel DBM compared with the predicate device, the in vivo implantation study summarized in this report was performed. Both the proposed membrane and the predicate membrane were evaluated for degradation rate after implantation over dental extraction sites which were filled with bone graft material.

Four (4) animals were euthanized at each time point: Day 7, Day 42, and Day 84. The time points were selected based on the degradation data available from previous in vivo experiments. The Day 7 animals were intended to provide the base line data, the Day 42 animals to provide data for intermediately degraded membranes and the Day 84 animals to provide data for (almost) completely degraded membranes.

This study was designed to answer the following questions raised during the review of a 510(k) premarket notification to FDA:

- Degradation Rate
Is there a comparable rate of resorption for both membranes? In case of slow resorption rates, does this increase the risk of late inflammation and infection?
- Tissue Integration
Does the implanted membrane cause unacceptable tissue reactions, both clinically (wound dehiscence or other complications) as well as histologically?
- Bone formation
Does tissue reaction to the membrane interfere with the bone formation process, i.e. influence the amount of bone formation at the adjacent implanted bone defect? The presence of comparable levels of bone formation in sites covered with the proposed and predicate devices will aid in the demonstration of substantial equivalence in terms of device performance.

5. ANIMAL MODEL

Dogs were used in this study because they have historically been used in safety evaluation studies and the guidelines have no other approved alternative (non-animal) methods. The animal species, number, and route of test material administration were selected taking into consideration the recommendations provided by Giannobile and Nevins 2011 (see References section).

The oral cavity of the canine most closely represents the human oral cavity in terms of potential for microbial load, dentition, and contact with the bone grafting material (Bio-Oss® spongiuous bone substitute).

The oral route of exposure was selected because it is the intended route of human exposure. The test or control article was placed over a bone graft (Bio-Oss® spongiuous bone substitute) after tooth extraction. The Bio-Oss bone substitute is a state of the art bone grafting material in GTR and GBR. In addition, in clinical practice the use of Bio-Oss is frequently combined with the Bio-Gide membrane.

Four premolar teeth, both mandibular and maxillary premolar 4s, were surgically extracted in their entirety, including both roots. This resulted in a bone defect of approximately 15-20 mm length by 5 mm width by 20 mm depth. Extraction sockets were cleaned of any remaining soft connective tissue. The gingiva surrounding the extraction sockets were undermined. Each extraction socket was filled with Bio-Oss bone substitute. The right maxilla and mandible extraction sockets were covered with the test article. The left maxilla and mandible extraction sockets were covered with the control article. The membranes were trimmed as necessary to allow an overlap of approximately 2 mm between the membrane and the defect walls. The membranes were inserted between the gingiva and the bone such that the filled extraction sockets were completely covered. The test article membrane was tacked to the gingiva with non-absorbable suture to prevent the membrane from dislodging in case of wound dehiscence. Each extraction socket was closed by closing the gingiva over the membrane with non-absorbable sutures.

\\SBS\Document Control\QP documents\QP-VER\QP-VER-03005\QP-VER-03005-A.docx

matricel	Verification Report	Document: QP-VER-03005-A
		Release: 27.Jan.14
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1	Page: 3 of 19

The animals were fed a soft diet post-operatively for approximately 7 days to prevent interference of daily nutrition related forces with the operation sites. Other protective measures are not practical in the canine model. The animals were transitioned to a solid diet between Day 7 and Day 14.

A detailed description of the experimental procedure is provided in Section 7 of the final GLP report.

6. EVALUATION CRITERIA

Radiographic Observation

Animals were pre-screened for the study with dental radiographs performed during the week prior to the initial teeth extraction and administration of the test and control article to ensure that the animals had adult dentition. Dental radiographs of the upper and lower left and right fourth premolars were also taken for all study animals post-operatively and on the day of necropsy. Test and control article postoperative and pre-necropsy dental radiographs were evaluated for obvious differences in signs of bone resorption and bone formation in the sockets of the mandible and maxilla irrespective of the time point.

Clinical Observation

After implantation, animals were observed daily for implant site healing (wound dehiscence or other complications) and for clinical signs of toxicity for the entire duration of the study. Findings were recorded as they were observed.

Histopathology

Histopathology was performed using hematoxylin/eosin and Masson's trichrome stained sections of the implant sites. Each lesion was listed and coded by the most specific topographic and morphologic diagnoses, severity, and distribution (see Table 1 below). A five step/severity grading system was used to define gradable lesions for comparison between treated and control groups. Representative images from the stained slides and examples of scoring are provided in ATTACHMENT B of the final GLP report.

Each transverse section was evaluated and graded separately. The following parameters were assessed histopathologically at the test article-graft interface, the test article-soft-tissue interface, the test article-bone interface, and at the graft-bone interface:

- foreign debris,
- degree of hemorrhage,
- intervening non-calcified material (test and control article),
- amount of adjacent new bone formation (immature/woven bone vs. mature/lamellar bone),
- resorption/remodeling (presence and activity of osteoclasts),
- polymorphic nuclear cells,
- mononuclear cells (lymphocytes/ macrophages),
- plasma cells,
- giant cells,
- vascularization,
- fibrosis,
- fatty infiltrate, and
- necrosis.

matricel	Verification Report	Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1	Release: 27.Jan.14
		Page: 4 of 19

Table 1: Details of the histopathological scoring system used.

Sites	Grading Scale	Terms
A= Upper left tooth, Rostral third of implant site.	0 = Normal	F = Focal
B= Upper left tooth, Middle region of implant site.	1 = Minimal	M = Multifocal
C= Upper left tooth, Caudal third of implant site.	2 = Mild	
D= Upper right tooth, Rostral third of implant site.	3 = Moderate	
E= Upper right tooth, Middle region of implant site.	4 = Marked	
F= Upper right tooth, Caudal third of implant site.	5 = Severe	
G= Lower left tooth, Rostral third of implant site.		
H= Lower left tooth, Middle region of implant site.		
I= Lower left tooth, Caudal third of implant site.		
J= Lower right tooth, Rostral third of implant site.		
K= Lower right tooth, Middle region of implant site.		
L= Lower right tooth, Caudal third of implant site.		

The upper and lower left sided sites (A-C and G-I) were control article implant sites. The upper and lower right sided sites (D-F and J-L) were test article implant sites. The term NA (not applicable) was used to denote when a particular section could not be evaluated due to lack of the proper tissue segments for evaluation. For example, when the Test or Control article was not present, the interface between that material and the bone graft material, soft tissue or bone could not be evaluated.

Histomorphometry

Histomorphometric analyses (membrane as well as bone and bone graft measurements) were performed on 3 Masson's trichrome stained sections from each site per dog. The histomorphometrist was blinded as to specimen identity when performing the measurements and analysis.

5. RESULTS

The most important results of the GLP study are summarized below.

Clinical Observations (Section 9.1 of the final GLP report)

Clinical observations were normal except for procedure related edema/swelling of the mandibular and maxillary regions for up to 6 days after the extractions. There were no findings of incisional dehiscence on clinical observations during the duration of the study, for neither test nor control membrane.

Summary of Radiographic Observations (Section 9.3 of the final GLP report)

All animals included in the study had adult dentition and were suitable for testing. There were no significant differences in the comparison of new bone formation between the left and right sides of dental radiographs as evident by comparing the radiographic opacity of the sockets of test and control articles. There was no evidence of bone resorption on any radiographs pre-necropsy for either test or control articles at any time points.

Histopathology (Section 9.4 of the final GLP report)

Note: In order to improve the readability of the histopathological results, the original scoring tables as provided in ATTACHMENT B of the final GLP report were re-organized to facilitate easier comparison between the data obtained for the test and control article. In addition, scorings that were not considered significant for all three necropsy times, due to the rare incidence (foreign debris, giant cells, fatty infiltrates and necrosis), were deleted from the re-organized tables. The re-organized tables are provided as attachments to this report.

	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 5 of 19

Day 7: No significant differences in the different aspects of the tissue reaction were noted between the control and test article groups at Day 7. Intervening non-calcified material (i.e. test or control article) had occasionally migrated away from the original implant site.

Day 42: No significant differences in the different aspects of the tissue reaction were noted between the control and test article groups at Day 42. The intervening non-calcified materials (test or control article) were at various stages of resorption, with smaller amounts of it noted at this time point compared to Day 7. In addition, the intervening non-calcified material had occasionally migrated away from the original implant site. Abundant new bone formation was noted at the areas of the extracted teeth. The bone grafts were undergoing active remodeling and were often attached to new woven bone.

Day 84: No significant differences in the different aspects of the tissue reaction were noted between the control and test article groups at Day 84. The intervening non-calcified materials (test or control article) were at advanced stages of resorption, with smaller amounts of it noted at this time point compared to Day 7. The amount and appearance of intervening non-calcified materials (test or control article) at Day 84 was similar to Day 42. In addition, the intervening non-calcified material (i.e. test or control article) had often migrated away from the original implant site. Abundant new bone formation was noted at the areas of the extracted teeth. The bone graft material was undergoing active remodeling and was often attached to new woven bone.

Histomorphology (Section 9.5 of the final GLP report):

Membrane Measurements: The individual data for membrane measurements at 7 days had values >0 at 7 days. At Day 42 and Day 84, there was only one membrane measurement (at 42 days) available, so no statistical analysis was performed on these data (see also the Discussion section below).

Bone and Bone Graft Measurements: In the membrane-adjacent ROI, due to the absence of membranes in almost all specimens collected after 7 days, there were no data to evaluate statistically except at 7 days, and there were no statistically significant differences between left and right at that time point.

In the membrane-distant ROI, Bone Area/Tissue Area increased with time and Bone Void Filler (BVF) Area/Tissue Area decreased with time (see Figure 1 below). There were no significant differences in total (mandibular + maxillary) values at any time point.

matricel	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 6 of 19


Figure 1: Bone Area/Tissue Area and Bone Void Filler (BVF) Area/Tissue Area for both test as well as control article, as a function of implantation time.

5. DISCUSSION

In this study, a canine extraction socket model was used to evaluate the test and control membrane for its degradation rate after implantation over dental extraction sockets. According to our knowledge, all reasonably possible measures were taken to prevent that normal animal behavior (eating, chewing, licking wounds, etc) would interfere with the stability of the implantation sites. However, these measures could not prevent that the implanted membranes had often migrated away from the original implant site, as described during histopathological evaluation.

In addition, artifacts caused by the histological preparation of the tissue sections were seen. Artifacts included splits or holes in the tissue sections, which appeared to represent tissue missing from the sections. In that case no bone histomorphometric measurements were made because it is impossible to determine how much tissue that should have been measured is missing. Pictures of every section are included in the Histomorphology Report (ATTACHMENT C in the final GLP report), so that this type of artifacts can be evaluated. As a consequence many measurements are missing from the histomorphometry, but these measurements are missing in a random fashion, so the remaining measurements are considered representative of the study results.

Migrated membranes and the artifacts caused by preparation of the tissue sections also explain the apparent difference in outcome of the membrane degradation data. Per the histomorphometry protocol, membranes and membrane-adjacent ROI were only measured if the membrane was present at the implantation site and can be reliably distinguished from the surrounding connective tissue. During histopathology, it was allowed to look at intervening non-calcified materials (test or control article) in the complete tissue section.

	Verification Report	Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1	Release: 27.Jan.14 Page: 7 of 19

During histomorphometry it was reported that tooth fragments were sometimes present in the tissue sections. These fragments resulted from the use of a bur during removal of any remaining root remnants and soft connective tissue from the extraction sockets. Although the extraction sockets were extensively rinsed after cleaning, this did not remove all the fragments.

7. CONCLUSIONS

This study was designed to answer the following questions raised during the review of a 510(k) premarket notification to FDA:

- Degradation Rate
Is there a comparable rate of resorption for both membranes? In case of slow resorption rates does this increase the risk of late inflammation and infection?
- Tissue Integration
Does the implanted membrane cause unacceptable tissue reactions, both clinically (wound dehiscence or other complications) as well as histologically?
- Bone Formation
Does tissue reaction to the membrane interfere with the bone formation process, i.e. influence the amount of bone formation at the adjacent implanted bone defect? The presence of comparable levels of bone formation in sites cover with the proposed and predicate devices will aid in the demonstration of substantial equivalence in terms of device performance.

The results of this study show that:

- Results - Degradation Rate
There were no differences in the degree of resorption or the extent of new bone formation between the test and control articles out to Day 84 following implantation.
- Results - Tissue Integration
There were no unacceptable tissue reactions observed during the study out to Day 84 following implantation. There were no differences in the inflammatory responses induced by the test and control articles implanted. There was no significant difference in the histological appearance between the test and the control articles after implantation. These results demonstrated that neither membrane induced any local toxic effect.
- Results - Bone Formation
There were no differences in the bone formation process when either the test article or the control article were used out to Day 84 following implantation

In summary, the results support that the Matricel DBM (test article) exhibits similar resorption and tissue reaction as compared to the Bio-Gide membrane (control article). On the basis of these results, for the purposes of the FDA 510(k) premarket notification, Matricel concludes that the Matricel DBM is substantially equivalent to the predicate Bio-Gide membrane.

ATTACHMENT 1													
Test/Control Article to Bone Graft Interface													
Day 7													
Jaw	Side	Article	Position	Animal/ Site	Hemorrhage	Intervening Non-Calified Material	Amount of New Bone Formation	Resorption/ Remodeling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	1D	1M	3F	0	0	0	0	0	1M	2M
				2D	2M	3F	0	0	3M	1M	0	2M	1M
				3D	2M	3F	0	0	0	1M	0	1M	1M
				4D	1M	3F	0	0	1M	0	0	1M	1M
			Middle	1E	2M	3F	0	0	0	0	0	1M	1M
				2E	1M	3F	0	0	0	1M	0	1M	1M
				3E	1M	3F	0	0	0	0	0	2M	1M
				4E	1M	3F	0	0	0	1M	0	2M	1M
			Caudal	1F	2M	3F	0	0	0	1M	0	1M	2M
				2F	2M	3F	0	0	0	1M	0	1M	1M
				3F	0	3F	0	0	0	1M	0	1M	1M
				4F	1M	3F	0	0	1M	1M	0	1M	1M
	Left	Control	Rostral	1A	3M	2F	0	0	3M	2M	0	0	1M
				2A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				3A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4A	2M	3F	0	0	1M	1M	0	2M	1M
			Middle	1B	2M	2F	0	0	3M	4M	0	0	1M
				2B	3M	3F	0	0	0	0	0	1M	2M
				3B	1M	3M	0	0	0	0	0	0	0
				4B	0	3F	0	0	1M	1M	0	1M	1M
			Caudal	1C	2M	3F	0	0	0	1M	0	0	0
				2C	4M	3F	0	0	0	0	0	0	0
				3C	1M	3F	0	0	1M	1M	0	2M	1M
				4C	1M	3F	0	0	1M	1M	0	1M	0
Lower	Right	Test	Rostral	1J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				3J	2M	3F	0	0	1M	1M	0	2M	1M
				4J	1F	1F	0	0	0	0	0	0	0
			Middle	1K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2K	1M	3F	0	0	0	0	0	0	0
				3K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4K	0	3F	0	0	1M	0	0	1M	1M
			Caudal	1L	0	2F	0	0	4M	1M	0	2M	1M
				2L	2M	3F	0	0	0	0	0	0	0
				3L	2M	3F	0	0	0	0	0	0	1M
				4L	2M	3F	0	0	3M	1M	0	1M	1M
	Left	Control	Rostral	1G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2G	1M	2F	0	0	0	0	0	0	1M
				3G	4M	3F	0	0	0	0	0	0	0
				4G	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	1H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2H	0	3F	0	0	2M	0	0	1M	1M
				3H	4M	3F	0	0	0	0	0	0	0
				4H	1M	3F	0	0	1M	1M	0	1M	1M
			Caudal	1I	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2I	1M	1M	0	0	0	0	0	0	0
				3I	2M	3F	0	0	2M	1M	0	2M	1M
				4I	1M	3F	0	0	0	0	0	1M	0

	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 9 of 19

ATTACHMENT 1													
Test/Control Article to Bone Graft Interface													
Day 42													
Jaw	Side	Article	Position	Animal/ Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/ Remodeling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	6D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13D	0	1F	3M	3M	0	0	0	0	2M
				14D	0	1F	2M	2M	0	0	0	1M	2M
			Middle	6E	0	1M	1M	0	0	0	0	1M	3M
				7E	0	1F	3M	3M	0	0	0	1M	2M
				13E	0	1M	4M	4M	0	0	0	1M	1M
				14E	0	1M	2M	1M	0	0	0	1M	2M
			Caudal	6F	0	1F	1M	1M	0	0	0	0	3M
				7F	0	1F	3M	3M	0	0	0	1M	2M
				13F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14F	0	1F	4M	3M	0	0	0	1M	1M
	Left	Control	Rostral	6A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7A	0	2F	0	0	1M	1M	0	1M	2M
				13A	0	1M	0	0	0	0	0	0	2M
				14A	0	2F	0	0	0	0	0	1M	2M
			Middle	6B	0	1F	0	0	0	1M	0	0	3M
				7B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13B	0	2F	3M	3M	0	0	0	1M	2M
				14B	0	2F	2M	2M	0	0	0	1M	2M
			Caudal	6C	0	1F	0	0	1M	0	0	1M	2M
				7C	0	2F	1M	2M	0	0	0	1M	3M
				13C	0	2M	4M	3M	0	1M	0	1M	1M
				14C	0	2F	2M	2M	0	0	0	1M	2M
Lower	Right	Test	Rostral	6J	0	1F	0	0	0	0	0	0	2M
				7J	0	1F	1F	0	0	0	0	2M	2M
				13J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14J	0	1F	1M	0	0	0	0	0	2M
			Middle	6K	0	2F	1F	0	0	0	0	0	2M
				7K	0	1F	0	0	0	1M	0	1M	2M
				13K	0	2M	1M	1M	0	1M	0	1M	2M
				14K	0	2F	2M	2M	0	0	0	1M	2M
			Caudal	6L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7L	0	1F	0	0	0	0	0	1M	2M
				13L	0	2F	0	0	0	1M	0	1M	2M
				14L	0	2F	1F	0	0	1M	0	1M	2M
	Left	Control	Rostral	6G	0	2F	2M	0	0	0	0	1M	2M
				7G	0	2F	0	0	0	0	0	1M	1M
				13G	0	2M	2M	2M	0	1M	0	0	1M
				14G	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	6H	0	2F	1M	1M	0	0	0	1M	2M
				7H	0	2F	2M	1M	0	0	0	1M	2M
				13H	0	2M	2M	1M	0	1M	1M	1M	1M
				14H	0	2F	0	0	0	1M	0	1M	1M
			Caudal	6I	0	2F	0	0	0	0	0	1M	2M
				7I	0	2F	0	0	0	0	0	2M	1M
				13I	0	2F	1M	1M	0	1M	1M	1M	1M
				14I	1M	1F	1F	0	0	0	0	1M	2M

matricel	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 10 of 19

ATTACHMENT 1													
Test/Control Article to Bone Graft Interface													
Day 84													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/Remodelling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	5D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				10D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12D	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5E	0	2F	2M	2M	0	0	0	1M	2M
				10E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12E	0	2F	0	0	0	0	0	1M	2M
			Caudal	5F	0	2M	2M	2M	0	1M	0	1M	2M
				10F	0	1M	2M	2M	0	1M	0	1M	2M
				11F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12F	0	2M	2M	1M	0	0	0	1M	2M
	Left	Control	Rostral	5A	0	1F	2M	3M	0	0	0	1M	2M
				10A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12A	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5B	0	2F	0	0	0	0	0	1M	2M
				10B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12B	0	1F	1F	1M	0	0	0	1M	2M
			Caudal	5C	0	2F	0	0	0	0	0	1M	2M
				10C	0	1F	1F	1M	0	1M	1M	1M	2M
				11C	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12C	0	1M	0	0	0	1M	0	1M	2M
Lower	Right	Test	Rostral	5J	0	1F	0	0	0	0	0	1M	2M
				10J	0	1M	0	0	0	0	0	1M	2M
				11J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12J	0	1F	1F	0	0	0	0	1M	2M
			Middle	5K	0	1F	0	0	0	0	0	1M	2M
				10K	0	1M	1M	1M	0	0	0	1M	2M
				11K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12K	0	1F	0	0	0	0	0	1M	2M
			Caudal	5L	0	1F	1F	1F	0	0	0	1M	1M
				10L	1M	2F	0	0	0	0	0	1M	2M
				11L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12L	0	1F	0	0	0	0	0	1M	2M
	Left	Control	Rostral	5G	1M	2F	0	0	0	0	0	1M	2M
				10G	0	1F	0	0	0	0	0	1M	2M
				11G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12G	0	1F	0	0	0	0	0	1M	2M
			Middle	5H	0	2F	0	0	0	0	0	1M	2M
				10H	0	1M	0	0	0	0	0	2M	2M
				11H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12H	1M	2F	1M	1M	0	0	0	1M	2M
			Caudal	5I	0	2F	0	0	0	0	0	1M	2M
				10I	0	2F	0	0	0	0	0	1M	2M
				11I	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12I	0	2F	0	0	0	0	0	1M	2M

ATTACHMENT 2													
Test/Control Article to Soft Tissue Interface													
Day 7													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calified Material	Amount of New Bone Formation	Resorption/Remodelling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	1D	2M	3F	0	0	0	1M	0	1M	1M
				2D	1M	2F	0	0	3M	1M	0	3M	1M
				3D	2M	3F	0	0	2M	1M	0	2M	1M
				4D	2M	3F	0	0	2M	1M	0	2M	1M
			Middle	1E	2M	3F	0	0	0	2M	0	2M	2M
				2E	2M	3F	0	0	1M	1M	0	2M	2M
				3E	1M	3F	0	0	1M	1M	0	1M	1M
				4E	1M	3F	0	0	1M	0	0	2M	1M
			Caudal	1F	1M	3F	0	0	2M	2M	0	3M	2M
				2F	2M	3F	0	0	1M	1M	0	1M	1M
				3F	1M	3F	0	0	0	0	0	0	0
				4F	2M	3F	0	0	2M	1M	0	2M	1M
	Left	Control	Rostral	1A	1M	2F	0	0	2M	1M	0	0	3M
				2A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				3A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4A	0	3F	0	0	2M	1M	0	1M	1M
			Middle	1B	3M	2F	0	0	3M	3M	0	3M	1M
				2B	1M	3F	0	0	0	1M	0	1M	1M
				3B	1M	3F	0	0	1M	1M	0	1M	1M
				4B	0	3F	0	0	1M	1M	0	1M	1M
			Caudal	1C	3M	3F	0	0	0	1M	0	2M	2M
				2C	3F	0	0	2M	2M	0	0	3M	1M
				3C	1M	3F	0	0	0	1M	0	1M	0
				4C	1M	3F	0	0	2M	1M	0	1M	1M
Lower	Right	Test	Rostral	1J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				3J	0	3F	0	0	0	0	0	0	0
				4J	1M	1F	0	0	2F	1F	0	0	0
			Middle	1K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2K	1M	3F	0	0	1M	0	0	1M	1M
				3K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4K	1M	3F	0	0	1M	0	0	1M	1M
			Caudal	1L	1M	2F	0	0	3M	1M	0	2M	1M
				2L	1M	3F	0	0	2M	1M	0	1M	1M
				3L	0	3F	0	0	0	0	0	0	1M
				4L	0	3F	0	0	2M	1M	0	2M	1M
	Left	Control	Rostral	1G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2G	1M	2F	0	0	1M	1M	0	0	1M
				3G	1M	3F	0	0	1M	1M	0	1M	0
				4G	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	1H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2H	1M	3F	0	0	1M	0	0	1M	1M
				3H	1M	3F	0	0	2M	1M	0	1M	1M
				4H	0	3F	0	0	1M	0	0	1M	M
			Caudal	1I	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2I	0	1M	0	0	2M	1M	0	1M	1M
				3I	1M	3F	0	0	1M	0	0	1M	1M
				4I	1M	3F	0	0	2M	0	0	1M	1M

ATTACHMENT 2													
Test/Control Article to Soft Tissue Interface													
Day 42													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/Remodelling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	6D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13D	0	2F	0	0	0	1M	1M	0	2M
				14D	0	1F	0	0	0	0	0	1M	2M
			Middle	6E	0	1M	0	0	0	0	0	0	2M
				7E	0	1F	0	0	0	0	0	1M	2M
				13E	0	1M	0	0	0	0	0	1M	2M
				14E	0	1M	0	0	0	1M	0	0	1M
			Caudal	6F	0	1F	0	0	0	1M	0	0	2M
				7F	0	1F	0	0	0	0	0	1M	3M
				13F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14F	0	1F	0	0	0	0	0	1M	2M
	Left	Control	Rostral	6A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7A	0	2F	0	0	0	0	0	1M	1M
				13A	0	1M	0	0	0	0	0	0	3M
				14A	0	3F	0	0	0	0	0	1M	2M
			Middle	6B	0	1F	0	0	0	1M	0	0	2M
				7B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13B	0	2F	0	0	0	0	0	1M	3M
				14B	0	2F	0	0	0	0	0	1M	2M
			Caudal	6C	0	1F	0	0	0	0	2M	1M	2M
				7C	0	2F	0	0	0	1M	1M	0	2M
				13C	0	2M	0	0	0	0	0	0	2M
				14C	0	2M	0	0	0	0	0	1M	2M
Lower	Right	Test	Rostral	6J	0	1F	0	0	0	0	0	0	2M
				7J	0	1F	0	0	1M	1M	0	2M	1M
				13J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14J	0	1F	0	0	0	0	0	1M	2M
			Middle	6K	0	2F	0	0	0	0	0	0	2M
				7K	0	1F	0	0	0	0	0	2M	2M
				13K	0	2M	0	0	0	0	0	1M	2M
				14K	0	2F	2M	2M	0	0	0	1M	2M
			Caudal	6L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7L	0	1F	0	0	0	0	0	1M	1M
				13L	0	2F	0	0	0	0	0	1M	2M
				14L	0	2F	1F	0	0	1M	0	1M	2M
	Left	Control	Rostral	6G	1M	1F	0	0	0	0	0	0	2M
				7G	0	2F	0	0	0	0	0	1M	2M
				13G	0	2M	0	0	0	1M	2M	1M	2M
				14G	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	6H	0	2F	0	0	1M	0	0	2M	2M
				7H	0	2F	0	0	0	0	0	0	2M
				13H	0	2M	0	0	0	0	0	1M	2M
				14H	0	2F	0	0	0	0	0	1M	1M
			Caudal	6I	0	2F	0	0	1M	1M	1M	1M	2M
				7I	0	2F	0	0	0	0	0	1M	2M
				13I	0	2F	0	0	0	1M	2M	1M	2M
				14I	0	1F	0	0	0	0	0	2M	2M

matricel	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 13 of 19

ATTACHMENT 2													
Test/Control Article to Soft Tissue Interface													
Day 84													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/Remodelling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	5D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				10D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12D	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5E	0	2F	0	0	0	0	0	1M	2M
				10E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12E	0	2F	0	0	0	0	0	1M	2M
			Caudal	5F	0	2M	0	0	0	0	0	2M	2M
				10F	0	1M	0	0	0	0	0	1M	2M
				11F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12F	0	2M	0	0	0	0	0	2M	2M
	Left	Control	Rostral	5A	0	1F	0	0	0	0	0	1M	2M
				10A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12A	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5B	0	2F	0	0	0	0	0	1M	2M
				10B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12B	0	1F	0	0	0	0	0	1M	2M
			Caudal	5C	0	2F	0	0	0	0	0	1M	2M
				10C	0	1F	0	0	0	1M	0	1M	2M
				11C	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12C	0	1M	0	0	0	0	0	1M	2M
Lower	Right	Test	Rostral	5J	0	1F	0	0	0	0	0	1M	2M
				10J	0	1M	0	0	0	0	0	0	2M
				11J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12J	0	1F	0	0	0	0	0	1M	2M
			Middle	5K	0	1F	0	0	0	0	0	1M	2M
				10K	1M	1M	0	0	0	0	0	1M	2M
				11K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12K	0	1F	0	0	0	0	0	1M	2M
			Caudal	5L	0	1F	0	0	0	0	0	1M	1M
				10L	0	2F	0	0	0	0	0	2M	2M
				11L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12L	0	1F	0	0	0	0	0	1M	2M
	Left	Control	Rostral	5G	2M	2F	0	0	1M	1M	1M	1M	2M
				10G	0	1F	0	0	0	0	0	1M	2M
				11G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12G	0	1F	0	0	0	0	0	1M	2M
			Middle	5H	0	2F	0	0	0	0	0	1M	1M
				10H	0	1M	0	0	0	0	0	1M	2M
				11H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12H	1M	2F	0	0	0	0	0	1M	2M
			Caudal	5I	0	2F	0	0	0	0	0	1M	2M
				10I	0	2F	0	0	0	0	0	2M	2M
				11I	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12I	0	2F	0	0	0	0	0	1M	2M

	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 14 of 19

ATTACHMENT 3													
Test/Control Article to Bone Interface													
Day 7													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/Remodeling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	1D	1M	3F	0	2M	0	1M	0	1M	2M
				2D	3M	2F	2M	3M	3M	1M	0	2M	1M
				3D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4D	1M	3F	0	1M	2M	1M	0	2M	1M
			Middle	1E	1M	3F	0	2F	0	1M	0	1M	1M
				2E	1M	3F	0	1M	0	1M	0	1M	1M
				3E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4E	0	3F	1M	3M	0	0	0	0	1M
			Caudal	1F	1M	3F	2M	2M	0	1M	0	1M	1M
				2F	2M	3F	0	3M	1M	1M	0	1M	1M
				3F	0	3F	0	3M	0	0	0	1M	0
				4F	1M	3F	2M	3M	1M	1M	0	2M	1M
	Left	Control	Rostral	1A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				3A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4A	0	3F	0	2M	1M	0	0	0	1M
			Middle	1B	2M	2F	0	3M	0	1M	0	1M	2F
				2B	0	3F	0	1M	0	1M	0	1M	1M
				3B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4B	0	3F	0	3M	0	1M	0	0	1M
			Caudal	1C	1M	3F	0	3M	0	2M	0	1M	3M
				2C	4M	3F	0	3M	0	1M	0	3M	1M
				3C	0	3F	1M	2M	1M	1M	0	2M	1M
				4C	0	3F	0	2M	0	0	0	0	1M
Lower	Right	Test	Rostral	1J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				3J	2M	3F	0	2M	0	1M	0	1M	1M
				4J	0	1F	0	2M	0	0	0	1M	1M
			Middle	1K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2K	2M	3F	1M	3M	0	0	0	1M	1M
				3K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				4K	0	3F	2M	1M	0	0	0	2M	1M
			Caudal	1L	0	2F	0	0	3M	1M	0	2M	2M
				2L	2M	3F	1M	1M	1M	1M	0	1M	1M
				3L	1M	3F	0	2M	1M	0	0	1M	1M
				4L	0	3F	2M	2M	1M	1M	0	1M	1M
	Left	Control	Rostral	1G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2G	0	2F	1M	1M	0	0	0	1M	1M
				3G	1M	3F	1M	2M	0	0	0	1M	2M
				4G	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	1H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2H	0	3F	1M	1M	1M	0	0	0	1M
				3H	1M	3F	2M	1M	0	1M	0	1M	1M
				4H	0	3F	0	3M	0	0	0	0	1M
			Caudal	1I	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2I	0	1M	0	1M	1M	0	0	2M	1M
				3I	0	3F	0	2M	0	1M	0	1M	1M
				4I	0	3F	0	1M	0	0	0	0	0

matricel	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 15 of 19

ATTACHMENT 3													
Test/Control Article to Bone Interface													
Day 42													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calified Material	Amount of New Bone Formation	Resorption/Remodeling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	6D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14D	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	6E	0	1F	0	2M	0	0	0	0	2M
				7E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14E	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Caudal	6F	0	1F	3M	3M	0	1M	0	0	1M
				7F	0	1F	2M	2M	0	0	0	1M	2M
				13F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14F	NA	NA	NA	NA	NA	NA	NA	NA	NA
	Left	Control	Rostral	6A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14A	0	3F	1M	2M	0	0	0	1M	2M
			Middle	6B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13B	0	2F	2M	2M	0	0	0	1M	2M
				14B	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Caudal	6C	0	1F	2M	2M	0	0	0	1M	2M
				7C	0	2F	3M	3M	0	0	0	0	2M
				13C	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14C	NA	NA	NA	NA	NA	NA	NA	NA	NA
Lower	Right	Test	Rostral	6J	0	1F	3M	3M	0	0	0	0	2M
				7J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14J	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	6K	0	2F	3M	4M	0	0	0	1M	1M
				7K	0	1F	0	3M	0	0	0	3M	2M
				13K	0	2M	0	0	0	0	0	1M	2M
				14K	0	2F	2M	2M	0	0	0	1M	2M
			Caudal	6L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7L	0	1F	2M	2M	0	0	0	1M	2M
				13L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14L	0	2F	1F	0	0	1M	0	1M	2M
	Left	Control	Rostral	6G	0	2F	2M	3M	0	0	0	3M	1M
				7G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				13G	0	2M	3M	2M	0	1M	0	1M	1M
				14G	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	6H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				7H	0	2F	0	2M	0	0	0	1M	2M
				13H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14H	0	1F	1M	2M	0	0	0	1M	2M
			Caudal	6I	0	2F	1M	2M	0	0	0	1M	1M
				7I	0	2F	0	1M	0	0	0	1M	1M
				13I	0	2F	0	2M	0	0	0	0	2M
				14I	0	1F	0	2M	0	0	0	1M	2M

	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 16 of 19

ATTACHMENT 3													
Test/Control Article to Bone Interface													
Day 84													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/Remodeling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	5D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				10D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11D	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12D	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				10E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12E	0	2F	1M	2M	0	0	0	1M	2M
			Caudal	5F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				10F	0	1M	2M	2M	0	0	0	1M	2M
				11F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12F	0	2M	2M	2M	0	0	0	1M	2M
	Left	Control	Rostral	5A	0	1F	2M	3M	0	0	0	0	2M
				10A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11A	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12A	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5B	0	2F	2M	3M	0	0	0	1M	2M
				10B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11B	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12B	0	1F	2M	2M	0	0	0	1M	2M
			Caudal	5C	0	2F	2M	3M	0	1M	0	1M	2M
				10C	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11C	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12C	0	1M	2M	2M	0	0	0	1M	2M
Lower	Right	Test	Rostral	5J	0	1F	2M	2M	0	0	0	1M	2M
				10J	0	1M	2M	3M	0	0	0	1M	2M
				11J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12J	0	1F	1M	1M	0	0	0	1M	2M
			Middle	5K	0	1F	2M	2M	0	0	0	1M	2M
				10K	0	1M	2M	2M	0	0	0	1M	2M
				11K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12K	0	1F	1M	1M	0	0	0	1M	2M
			Caudal	5L	0	1F	1M	1M	0	0	0	1M	1M
				10L	0	2F	3M	3M	0	0	0	0	2M
				11L	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12L	0	1F	2M	3M	0	0	0	1M	2M
	Left	Control	Rostral	5G	0	2F	0	0	0	0	0	1M	2M
				10G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12G	0	0	1M	1M	0	0	0	1M	2M
			Middle	5H	0	2F	1M	1M	0	0	0	1M	2M
				10H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11H	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12H	0	0	1M	1M	0	0	0	0	1M
			Caudal	5I	0	2F	1M	2M	0	0	0	1M	2M
				10I	0	2F	2M	2M	0	0	0	1M	2M
				11I	NA	NA	NA	NA	NA	NA	NA	NA	NA
				12I	0	0	1M	1M	0	0	0	1M	2M

matricel	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 17 of 19

ATTACHMENT 4													
Bone Graft to Bone Interface													
Day 7													
Jaw	Side	Article	Position	Animal Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/Remodelling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibroblasts
Upper	Right	Test	Rostral	1D	1M	0	0	2M	0	0	0	1M	2M
				2D	0	0	2M	3M	0	0	0	2M	1M
				3D	0	0	3M	2M	0	0	0	1M	2M
				4D	0	0	2M	2M	0	0	0	3M	2M
			Middle	1E	2M	0	2M	3M	0	0	0	0	2M
				2E	0	0	2M	2M	0	0	0	1M	1M
				3E	0	0	3M	2M	0	0	0	4M	2M
				4E	0	0	1M	3M	0	0	0	1M	1M
			Caudal	1F	0	0	3M	2M	0	0	0	1M	1M
				2F	0	0	2M	3M	0	0	0	1M	1M
				3F	0	0	3M	4M	0	0	0	1M	1M
				4F	0	0	3M	4M	0	0	0	1M	1M
	Left	Control	Rostral	1A	0	0	3M	2M	0	0	0	0	3M
				2A	2M	0	2M	1M	0	0	0	1M	2M
				3A	0	0	3M	3M	0	0	0	1M	1M
				4A	0	0	3M	3M	0	0	0	2M	1M
			Middle	1B	2M	0	2M	3M	0	0	0	2M	3M
				2B	2M	0	2M	2M	0	0	0	1M	2M
				3B	0	0	3M	3M	0	0	0	1M	2M
				4B	0	0	3M	3M	0	0	0	0	1M
			Caudal	1C	0	0	3M	3M	0	0	0	1M	2M
				2C	0	0	2M	3M	0	0	0	1M	1M
				3C	0	0	2M	2M	0	0	0	0	1M
				4C	2M	0	2M	2M	0	0	0	1M	2M
Lower	Right	Test	Rostral	1J	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2J	0	0	0	1M	0	0	0	1M	1M
				3J	0	0	1M	1M	0	0	0	1M	2M
				4J	3M	0	2M	3M	0	0	0	3M	2M
			Middle	1K	NA	NA	NA	NA	NA	NA	NA	NA	NA
				2K	0	0	3M	3M	0	0	0	1M	1M
				3K	0	0	1M	2M	0	0	0	1M	1M
				4K	0	0	2M	1M	0	0	0	1M	1M
			Caudal	1L	1M	0	2M	2M	0	0	0	0	2M
				2L	0	0	2M	1M	0	0	0	1M	1M
				3L	0	0	1M	2M	0	0	0	1M	1M
				4L	3M	0	2M	2M	0	0	0	1M	2M
	Left	Control	Rostral	1G	1M	0	2M	3M	0	0	0	2M	1M
				2G	0	0	1M	1M	0	0	0	1M	1M
				3G	0	0	3M	2M	0	0	0	1M	1M
				4G	0	0	2M	3M	0	0	0	1M	1M
			Middle	1H	0	0	3M	3M	0	0	0	1M	1M
				2H	0	0	2M	2M	0	0	0	1M	1M
				3H	0	0	3M	2M	0	0	0	1M	1M
				4H	0	0	3M	3M	0	0	0	1M	1M
			Caudal	1I	1M	0	2M	2M	0	0	0	2M	1M
				2I	0	0	3M	2M	0	0	0	1M	0
				3I	0	0	1M	2M	0	0	0	1M	2M
				4I	0	0	2M	3M	0	0	0	1M	2M

	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 18 of 19

ATTACHMENT 4													
Bone Graft to Bone Interface													
Day 42													
Jaw	Side	Article	Position	Animal/ Site	Hemorrhage	Intervening Non-Calcified Material	Amount of New Bone Formation	Resorption/ Remodelling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	6D	0	0	1M	3M	0	0	0	0	3M
				7D	0	0	4M	3M	0	0	0	1M	1M
				13D	0	0	4M	4M	0	0	0	1M	2M
				14D	0	0	4M	3M	0	0	0	2M	1M
			Middle	6E	0	0	3M	4M	0	0	0	1M	1M
				7E	0	0	4M	3M	0	0	0	1M	2M
				13E	0	0	4M	4M	0	0	0	2M	1M
				14E	0	0	4M	4M	0	0	0	1M	1M
			Caudal	6F	0	0	3M	4M	0	0	0	1M	1M
				7F	0	0	4M	3M	0	0	0	1M	1M
				13F	NA	NA	NA	NA	NA	NA	NA	NA	NA
				14F	0	0	4M	3M	0	0	0	1M	1M
	Left	Control	Rostral	6A	0	0	4M	4M	0	0	0	1M	4M
				7A	0	0	4M	3M	0	0	0	1M	2M
				13A	0	0	3M	3M	0	0	0	1M	1M
				14A	0	0	3M	4M	0	0	0	2M	1M
			Middle	6B	0	0	4M	3M	0	0	0	3M	2M
				7B	0	0	4M	4M	0	0	0	1M	1M
				13B	0	0	4M	3M	0	0	0	2M	1M
				14B	0	0	2M	3M	0	0	0	1M	2M
			Caudal	6C	0	0	1M	2M	2M	1M	1M	1M	2M
				7C	0	0	4M	3M	0	0	0	1M	2M
				13C	0	0	4M	3M	0	0	0	2M	2M
				14C	0	0	4M	3M	0	0	0	2M	1M
Lower	Right	Test	Rostral	6J	0	0	3M	4M	0	0	0	1M	1M
				7J	0	0	4M	4M	0	0	0	2M	2M
				13J	0	0	3M	2M	0	0	0	1M	1M
				14J	0	0	3M	4M	0	0	0	2M	2M
			Middle	6K	0	0	4M	4M	0	0	0	1M	1M
				7K	0	0	3M	3M	0	0	0	2M	1M
				13K	0	0	3M	3M	0	0	0	1M	1M
				14K	0	2F	2M	2M	0	0	0	1M	2M
			Caudal	6L	0	0	4M	4M	0	0	0	1M	2M
				7L	0	0	2M	2M	0	0	0	1M	1M
				13L	0	0	3M	3M	0	0	0	1M	2M
				14L	0	2F	1F	0	0	1M	0	1M	2M
	Left	Control	Rostral	6G	0	0	3M	4M	0	0	0	1M	1M
				7G	0	0	4M	4M	0	0	0	1M	2M
				13G	0	0	4M	2M	0	0	0	2M	1M
				14G	0	0	4M	4M	0	0	0	1M	1M
			Middle	6H	0	0	3M	4M	0	0	0	1M	1M
				7H	0	0	3M	1M	0	0	0	1M	2M
				13H	0	0	3M	2M	0	0	0	1M	1M
				14H	0	0	2M	3M	0	0	0	2M	1M
			Caudal	6I	0	0	2M	4M	0	0	0	1M	1M
				7I	0	0	1M	4M	0	0	0	2M	1M
				13I	0	0	3M	2M	0	0	0	1M	1M
				14I	0	0	4M	3M	0	0	0	2M	1M

	Verification Report		Document: QP-VER-03005-A
	Dental Barrier Membrane Implantation in Beagle Dogs Summary of GLP study 13-02011-G1		Release: 27.Jan.14
			Page: 19 of 19

ATTACHMENT 4													
Bone Graft to Bone Interface													
Day 84													
Jaw	Side	Article	Position	Animal/ Site	Hemorrhage	Intervening Non-Calified Material	Amount of New Bone Formation	Resorption/ Remodeling of Bone	Polymorphonuclear Cells	Mononuclear Cells	Plasma Cells	Vascularization	Fibrosis
Upper	Right	Test	Rostral	5D	0	0	3M	3M	0	0	0	1M	2M
				10D	0	0	2M	2M	0	0	0	1M	2M
				11D	0	0	1M	2M	0	0	0	1M	3M
				12D	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Middle	5E	0	0	3M	3M	0	0	0	2M	2M
				10E	NA	NA	NA	NA	NA	NA	NA	NA	NA
				11E	0	0	2M	2M	0	0	0	2M	3M
				12E	0	0	1M	1M	0	0	0	1M	2M
			Caudal	5F	0	0	2M	2M	0	0	0	2M	1M
				10F	0	0	3M	2M	0	0	0	3M	2M
				11F	0	0	2M	3M	0	0	0	1M	2M
				12F	0	0	1M	1M	0	0	0	2M	2M
	Left	Control	Rostral	5A	0	0	3M	4M	0	0	0	2M	2M
				10A	0	0	3M	4M	0	0	0	1M	2M
				11A	0	0	3M	3M	0	0	0	2M	2M
				12A	0	0	3M	2M	0	0	0	1M	3M
			Middle	5B	0	0	3M	3M	0	0	0	1M	2M
				10B	0	0	3M	4M	0	0	0	2M	3M
				11B	0	0	2M	3M	0	0	0	2M	2M
				12B	0	0	4M	3M	0	0	0	3M	2M
			Caudal	5C	0	0	4M	3M	0	0	0	2M	3M
				10C	0	0	2M	2M	0	0	0	1M	2M
				11C	0	0	2M	2M	0	0	0	2M	3M
				12C	0	0	2M	2M	0	0	0	2M	2M
Lower	Right	Test	Rostral	5J	0	0	3M	2M	0	0	0	1M	2M
				10J	0	0	3M	3M	0	0	0	2M	1M
				11J	0	0	2M	2M	0	0	0	1M	1M
				12J	0	0	3M	2M	0	0	0	2M	1M
			Middle	5K	0	1F	1M	3M	0	0	0	2M	2M
				10K	0	0	4M	4M	0	0	0	2M	3M
				11K	0	0	2M	2M	0	0	0	1M	1M
				12K	0	0	2M	2M	0	0	0	1M	2M
			Caudal	5L	0	0	3M	4M	0	0	0	2M	2M
				10L	0	0	4M	4M	0	0	0	2M	3M
				11L	0	0	3M	3M	0	0	0	2M	2M
				12L	0	0	3M	3M	0	0	0	2M	1M
	Left	Control	Rostral	5G	NA	NA	NA	NA	NA	NA	NA	NA	NA
				10G	0	0	2M	2M	0	0	0	1M	2M
				11G	0	0	3M	3M	0	0	0	2M	2M
				12G	0	0	2M	2M	0	0	0	2M	2M
			Middle	5H	0	0	2M	2M	0	0	0	2M	1M
				10H	0	0	4M	4M	0	0	0	3M	2M
				11H	0	0	3M	3M	0	0	0	1M	3M
				12H	NA	NA	NA	NA	NA	NA	NA	NA	NA
			Caudal	5I	0	0	3M	3M	0	0	0	1M	2M
				10I	0	0	5M	4M	0	0	0	2M	3M
				11I	0	0	2M	2M	0	0	0	1M	2M
				12I	0	0	3M	3M	0	0	0	1M	2M

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