Operational Documentation

Helitene® ABSORBABLE COLLAGEN HEMOSTATIC AGENT

DESCRIPTION

HELITENE® Absorbable Collagen Hemostatic Agent is a soft, white, pliable absorbent hemostatic agent, supplied in a fibrillary form.

The basic material from which the HELITENE hemostatic agent is fabricated is collagen obtained from bovine deep flexor (Achilles) tendon. The tendon is known to be one of the purest sources of collagen that can be readily obtained and processed in commercial amounts. The HELITENE hemostatic agent being derived from this tendon, is expected to be very consistent material. Because of the initial purity of the collagen source and the further purification steps during processing of the HELITENE hemostatic agent, the practitioner can expect uniform behavior from this topical hemostat from one application to the next.

INDICATIONS

The HELITENE hemostatic agent is indicated in surgical procedures as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

POTENTIAL CONSUMERS AND OPERATING CONDITIONS

The products must be used by qualified personnel in the operating room

INFORMATION FOR USAGE

On contact with blood, collagen is known to cause aggregation of platelets. Platelets deposit in large numbers on the collagen structure, degranulate, and release coagulation factors that, together with plasma factors, enable the formation of fibrin. The structure of the HELITENE hemostatic agent provides a three dimensional matrix for the additional strengthening of the clot.

The HELITENE hemostatic agent effectively controls bleeding usually within two to five minutes when applied directly to the bleeding site. Excess the HELITENE hemostatic agent should be removed from the site after hemostasis is achieved.

The HELITENE hemostatic agent Absorbable Collagen Hemostatic Agent is designed to be totally absorbable if left in situ after hemostasis. If desired, the HELITENE hemostatic agent may be recovered after hemostasis is accomplished using dry forceps or irrigation. Implant studies in animals have demonstrated the HELITENE Absorbable Collagen Hemostatic Agent to be absorbed with tissue reaction similar to that observed with other absorbable hemostatic agents.

The collagen hemostatic agent absorption was evaluated after subcutaneous and intrahepatic implantation in rats. In one out of live animals, complete subcutaneous absorption was observed by day 14, and by day 56, three out of four animals had complete absorption. Complete intraperitoneal absorption was not observed by day 56.

As shown with other hemostatic agents, the implantation of the HELITENE hemostatic agent also elicits a similar foreign body reaction.

The HELITENE Absorbable Collagen Hemostatic Agent has been evaluated in vitro for the enhancement of bacterial growth of Staphylococcus aureus and Escherichia coli. Enhancement of bacterial growth did not occur for either organism.

In vivo studies using guinea pigs showed that incidence of infection (abscess) of incision sites inoculated with Staphylococcus aureus was not enhanced by the presence of the collagen hemostatic agent when

compared to another collagen hemostatic agent. However, extent of wound infection tended to be greater than control with the HELITENE Absorbable Collagen Hemostatic Agent and another collagen hemostatic agent tested. This tendency is observed with many foreign substances.

The HELITENE Absorbable Collagen Hemostatic Agent was evaluated for potential allergenic sensitivity. A guinea pig maximization study showed that it did not produce irritation or contact sensitization. A chemical assay of the HELITENE Absorbable Collagen Hemostatic Agent compared to one other collagen hemostat showed significantly less specific glycoprotein immunoreactive substances. A hemagglutination study was conducted evaluating the HELITENE Absorbable Collagen Hemostatic Agent as the antigen. There was no agglutination observed.

CONTRAINDICATIONS

The HELITENE hemostatic agent should not be used for patients with a known history of hypersensitivity to bovine derived products.

PRECAUTIONS

It has been shown that fragments of microfibrillar hemostatic agents may be passed through the filter of blood saving devices and returned to the patient. The use of microfibrillar collagen in patients undergoing surgical procedures involving blood saving devices should be avoided.

As with other hemostatic agents, it is not recommended that the HELITENE hemostatic agent be left in an infected or contaminated space.

The HELITENE hemostatic agent is not intended to be used to treat systemic coagulation disorders.

Only the amount of the HELITENE hemostatic agent necessary to produce hemostasis should be used. After approximately 10-15 minutes, excess material should be removed. This is usually possible by removing the HELITENE hemostatic agent using dry forceps or irrigation. In otolaryngological surgery, precaution against aspiration should include removal of excess dry material.

Clinical studies were not conducted in ophthalmological or urological surgical procedures, therefore, the risk to benefit ratio should be assessed prior to use in these procedures.

There are no well controlled studies in pregnant women; therefore, the HELITENE hemostatic agent should be used in pregnant women only when the benefit outweighs the risk.

SINGLE-USE DEVICE

The HELITENE hemostatic agent is supplied in a single-use package and is guaranteed to be sterile and non-pyrogenic unless opened or damaged. The product is intended for use as an absorbable implant and is not to be reused. Reuse of the device can result in contamination and/or disease transmission. Any attempt to resterilize or reuse the product/components will damage the matrix and impair its ability to function as intended. All unused pieces must be discarded.

ADVERSE REACTIONS

Adverse reactions reported with another microfibrillar collagen agent that were possibly related to its use were adhesion formation, allergic reaction, foreign body reaction, and subgaleal seroma (report of a single case). The use of microfibrillar collagen in dental extraction sockets has been reported to increase the incidence of alveolgia.

Other microfibrillar collagen have been reported to cause interference with the healing of skin edges when used in the closure of skin incisions and to reduce the strength of methyl-methacrylate adhesive when used to attach prosthetic devices to bone surfaces. Transient laryngospasm due to aspiration of dry material has been reported following the use of another microfibrillar collagen in tonsillectomy procedures.

Since the HELITENE hemostatic agent is a collagen-based product, adverse reactions experienced with other collagen hemostatic agents may be related.

ADMINISTRATION

The HELITENE hemostatic agent should be placed directly on the bleeding surface with pressure applied using a dry gauze. The period of time necessary to apply pressure will vary with the type and amount of bleeding to be controlled. In general, one to five minutes should be sufficient. It has been shown that hemostasis usually occurs within two to five minutes. The amount of the HELITENE hemostatic agent necessary to achieve hemostasis will depend on the nature and amount of bleeding to be controlled. In surgical procedures the pad can be used by applying small squares to bleeding areas. Either side of the HELITENE hemostatic agent pad may be applied to the bleeding site.

The HELITENE hemostatic agent will adhere to wet gloves, instruments and tissue surfaces. Dry forceps should be used to apply the HELITENE hemostatic agent to facilitate handling and placement of the material.

The HELITENE hemostatic agent may be left in situ whenever necessary. However, the surgeon should remove any excess the HELITENE hemostatic agent prior to wound closure.

Possibility of removal: Prior to wound closure. After About 10-15 minutes, the excess material should be removed. Usually hemostatic HELITENE can be removed by using dry forceps or irrigation.

In otolaryngological surgery, precautions should be taken to avoid aspiration when removing excess dry material.

Interaction of implant materials with other materials and substances: The manufacturer does not provide the possibility of interaction with drugs/materials/substances

Information about the methods of determining the position of the implant: Not applicable.

INFORMATION ABOUT THE AVAILABILITY OF MEDICAL DEVICE MATERIALS OF ANIMAL ORIGIN

Integra LifeSciences Corporation manufactures Helitene® ABSORBABLE COLLAGEN HEMOSTATIC AGENT from collagen – Integra Collagen. Integra Collagen is produced by Integra LifeSciences Corporation from bovine deep flexor tendon sourced from the United States and New Zealand. Each shipment of tendon has I Certificate of Conformance. Bovine tendon is subjected to enzymatic and alkali treatments, resulting in highly purified Type I collagen – Integra Collagen – that is non-viable and has been inactivated for viruses and TSEs.


SYMBOLS USED ON LABELING

	Consult Instructions for Use
	Expiration date
	Do not re-use
LOT	Lot number
	Do Not Use if package is damages
EC REP	Authorized Representative in the European Community
STERILE EO	Sterilized Using Ethylene Oxide
	Product compiles with requirements of directive 93/42/EEC
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or the order of a physician or practitioner
	Manufacturer
REF	Catalog number
	Do Not re-sterilize
	Peel here

Non-pyrogenic - indicating that the medical device non-pyrogenic

Integra LifeSciences Corpotation
311 Enterprise Drive
Playnsboro, NJ 08536
Telephone: 800-654-2873
Outside USA: 609-275-0500
Fax: 609-275-5363
www.integralife.com

Integra LifeSciences Services (France) SAS
Immeuble Sequoia 2
97 allée Alexandre Borodine
Parc Technologique de la Porte des Alpes
69800 Saint Priest, France
Telephone: 33 (0) 4 37 47 59 10
Fax: 33 (0) 4 37 47 59 29

HOW SUPPLIED

Sterile HELITENE Absorbable Collagen Hemostatic Agent-Fibrillar Form is supplied in the following sizes:

REF	Size	Quantity
1692-ZLI	0.5 gram	5/box
1692-ZKI	1.0 gram	5/box

Contents of the package are guaranteed sterile and non-pyrogenic unless the package is opened or damaged. Avoid excessive heat and humidity.

Storage and transportation conditions

	Environmental conditions (storage)	Environmental conditions (during transportation)
Temperature range	15-30 °C	15-30 °C

Relative humidity	No more than 65%	No more than 65%
-------------------	------------------	------------------

Special time intervals for transportation are not provided

UTILIZATION

Materials must be disposed of and/or recycled in accordance with applicable laws.

Waste class:

for products with expired class - A;

for products used class - B.

Local Representative in Russia

Neuroproject Ltd

127486, Moscow, Russia, Korovinskoye shosse, 10, bld. 2, floor 1, room 101

Tel. +7 495 6608173


e-mail: info@neurop.ru

Authorized Signature:

Name: Diana Bordon

Title: Director Regulatory Affairs

INTEGRA

Signature: 

Date: 03 Oct 2019

Jennifer P. Daniels
NOTARY PUBLIC OF NEW JERSEY
ID# 2285997, My Commission Expires 4/1/2021

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