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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SCENKER BIOLOGICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

ShaoYunjie

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(Date: Sep. 29, 2019)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY

SCENKER Biological Technology
Co., Ltd.

CUI GE/Legal

representative

19 September, 2019

Culture medium bottles IOHOHa[®]

for automated bacteriological systems IOHOHa[®] LABSTAR
in models

INSTRUCTION FOR USE

Foreword

We thank you for purchasing Culture medium bottles for LABSTAR Automated bacteriological Systems.

Bottles are designated for cultivation and detection of pathogens from clinical blood and other normally sterile human body fluids. The bottles are used with LABSTAR **ЮХОНА**[®] Automated bacteriological Systems.

The Bottles are ready for use, and are designated for single-use.

It is recommended to use the Bottles and the System in laboratory which corresponds to the second level of biological safety. Measures of biological protection shall be assumed in course of its operation. Used bottles shall be sterilized by means of autoclaving or other safe method, provided by local legislation requirements.

Product data:

Permit to manufacture No.: L.S.Y.J.X.S.C.X. No. 20110052

Manufacture date: specified on the label.

Product standard: YZB/ LU 0272-2013 dated 01 July 2013.

First instruction issue date: 01 July 2013.

Instruction version for registration in the Russian Federation:
v.2 dated 19 September 2019.

The product will be presented under **ЮХОНА**[®] trademark on the Russian Federation market.

Scenker Biological Technology Co., Ltd.

Index

1. Information on medical product.....	4
1.1. Medical product name.....	4
1.2. Information on medical product manufacturer	4
1.3. Information on medical product authorized representative	4
2. Description and operation	4
2.1. Brief information on medical product	4
2.2. Medical device technical specifications and characteristics, guaranteed by manufacturer	5
2.3. Composition (delivery completeness) of the medical product	15
2.4. Warnings, safety measures and restrictions during medical product use	15
2.5. Package, marking and warning designations description	16
3. Medical product operation	22
3.1. Adjustment, assemblage, setting up, gauging and other activities, required to put the product into operation.....	22
3.2. Operation procedure.....	23
3.3. Quality control	28
3.4. Medical product application risks, counterindications, expected and predictable side effects, connected to medical product use as intended.....	30
4. Information on maintenance and current repairs	30
5. Storage and transportation.....	30
5.1. Transportation	30
5.2. Storage.....	30
6. Warranty liabilities	31
7. Medical product disposal and liquidation procedure.	31
Appendix A	33
Appendix B	34
Appendix C	37

1. Information on medical product

1.1. Medical product name

Culture medium bottles **ЮHOHA®** for automated bacteriological systems
ЮHOHA® LABSTAR in models (hereinafter referred as medium bottle(s)):

1. Culture medium Bottle **ЮHOHA®** for aerobe cultivation;
2. Antibiotic-neutralized culture medium bottle **ЮHOHA®** for aerobe cultivation;
3. Antibiotic-neutralized culture medium bottle **ЮHOHA®** for anaerobe cultivation;
4. Antibiotic-neutralized culture medium bottle **ЮHOHA®** for aerobe cultivation for children.

1.2. Information on medical product manufacturer

SCENKER Biological Technology Co., Ltd., East End of Weier Road,
Fenghuang Industrial Park, Liaocheng, Shandong Province, China, 252000.

1.3. Information on medical product authorized representative

You shall refer to the manufacturer's authorized representative in the Russian Federation regarding quality and maintenance issues of medical products "Automated bacteriological System **ЮHOHA®** LABSTAR, design versions: **ЮHOHA®** LABSTAR 50, **ЮHOHA®** LABSTAR 100 with the accessories" and "Culture medium bottles **ЮHOHA®** for automated bacteriological systems **ЮHOHA®** LABSTAR in models", manufactured by SCENKER Biological Technology Co., Ltd, China:

Triopromtrade, Limited Liability Company (Triopromtrade, LLC),
The Russian Federation, 620028, Sverdlovsk region, Yekaterinburg, 28
Kirova st./43 Toledova st. Tel.: +7 (343) 270-75-00, Fax: +7 (343) 246-91-56.

E-mail: info.triopromtrade@mail.ru

2. Description and operation

2.1. Brief information on medical product

2.1.1. Medical product designation, provided by manufacturer

Culture medium bottles are designated for cultivation and detection of pathogens from clinical blood and other normally sterile human body fluids. The bottles are used with automated bacteriological systems **ЮHOHA®** LABSTAR.

2.1.2. Medical product applicable scope and potential consumers

This is used in medical institutions, clinical diagnostic centers. Users are qualified medical technicians, occupied in microbiological researches conduction.

2.1.3. Risk class and applicable classification rules

Medical product usage potential risk class in accordance with medical products classification nomenclature is 2a (i. 9.4 of the Order of Ministry of Health of Russia dated 06.06.2012 No.4H).

According to the Guideline “On medical products”, issued in 2002 by the Chinese Department for Medical Products and Drugs Control, the product is classified as “Reagents for laboratory in vitro diagnostics 6840” and belongs to the second class of medical products.

2.1.4. Indications for use

The necessity of cultivation and detection of pathogenic microorganisms in blood samples and other sterile biological fluids of the human body during microbiological studies to support the diagnosis of infectious diseases

2.1.5. Counterindications

The bottles are used only for in vitro diagnostics.

2.1.6. Possible side effects

Not applicable

2.2. Medical device technical specifications and characteristics, guaranteed by manufacturer

2.2.1. Models

“Culture medium bottles **ЮHOHA®** for automated bacteriological systems **ЮHOHA®** LABSTAR in models” have four models depending on pathogens type and clinical examinations scope:

- Culture Medium Bottle **ЮHOHA®** for aerobe cultivation (color of label and removable cover is red, designation ) – a standard culture medium bottles designed for adult patients not treated with antibiotics;

- Antibiotic-neutralized Culture Medium Bottle **ЮОНОН**[®] for aerobe cultivation (color of label and a removable cover is blue, designation ■), culture medium bottles designed for adult patients treated with antibiotics;
- Antibiotic-neutralized Culture Medium Bottle **ЮОНОН**[®] for anaerobe cultivation (color of label and removable cover is cadmium orange, designation □), culture medium bottles designed for adult patients treated with antibiotics;
- Antibiotic-neutralized Culture Medium Bottle **ЮОНОН**[®] for aerobe cultivation for children (the color of the label and removable cover is pink, designation ●), culture medium bottles designed for children.

Culture medium bottles are compatible with all models of automated bacteriological systems **ЮОНОН**[®] LABSTAR.

Image 1. Models

2.2.2. Description

Once the sample culture has been completed, culture medium bottles shall be put in the System for cultivation and continuous monitoring whether there is any growth of microorganisms in the Bottle by means of fluorescent method.

The optical sensitive element, located on the bottles bottom (polymer sensor), which starts shining under the influence of substance of metabolism that is generated by bacteria. If there are microorganisms in the sample, which was put in the Bottle, they produce carbon dioxide in course of substrates growth and disposal. LED gives off light, and this light reflects from the sensor, photoelectric detector identifies such light, and sends signal to the system computer analytic block, where the analogue signal is transformed into the discrete one. Due to

carbon dioxide accumulation the sensor reflects more light. Cвет, The results shall be compared to initial sensor indexes. If CO₂ content is high, the sample shall be considered to be positive, if CO₂ level doesn't change after incubation in course of the time specified by operator, the sample shall be considered to be negative. The visual reporting system activates on the monitor once positive samples have been detected.

Based on the measurement data system computer analytical block creates the growth curve, which view allows to draw a conclusion on pathogens presence or lack in this sample. The rate of growth of pathogens is directly proportional to the volume of the biomaterial (in the range recommended for each type of bottles), the optimal detection of pathogens is achieved with the maximum volume of the sample. On the diagram you can see the point of triggering an alarm in case of the presence of the positive sample.


Image 2 Growth curve of the positive sample


Image 3 Growth curve of the negative sample

The bottles are ready for use, the bottles shall be used once.

2.2.3. List of detected microorganisms

Culture medium bottles are designated for cultivation and detection of the following microorganisms in blood sample and other normally sterile human body fluids:

- Aerobic bacteria:
- Gram (+) Coccus:
 - Streptococcus pneumoniae
 - Streptococcus pyogenes
 - Staphylococcus aureus
 - Alcaligenes faecalis
- Gram (-) Coccus:

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- Fungus

- *Candida albicans*

The list of diseases, caused by these microorganisms, is introduced in Appendix B.

2.2.4. Information on materials

Culture medium bottles are packed leak-proof with the cork. The Bottle's body material is polyethylene terephthalate resin. There is a coloured label on the body.


Image 4. Culture medium bottle

1 – mixed cap with removable part, 2 – cork, 3 – bar code, 4 – filling marks, 5 – culture medium, 6 – resin granules, 7 – polymer sensor

Combined cap with removable part (1) – guarantees the Bottle's cleanness (but not sterility), it is manufactured from the aluminum stripe and removable part made of coloured plastic of the same colour the label is, which simplifies identification.

Cork (2) – butyl cork is fixed by the aluminum cap.

Bar code (3) – there is a unique bar code on each Bottle label. Bar code detachable section shall be stuck in the laboratory register or other document. Bar code information shall be scanned before placing Bottles in the System. Putting bar code in the System is required for the correct algorithm of processing of data of measurements in the Bottle, and also for putting the Bottle in a corresponding mesh for the following cultivation.

Filling marks (4) – 5 ml marks are designated for approximate determination of the volume of sample, put in the Bottle.

Culture medium (5) – yellowish liquid with granules of adsorption resin* intended for adsorption and neutralization of antibiotics (6). Culture medium volume in the Bottle amounts to 25 ml.

* exclusive model Culture Medium Bottle ЮНННН[®] for aerobe cultivation.

Polymer sensor (7) – affected by the bacteria substance of metabolism it starts to shine, the strength of shining increases as the quantity of bacteria increases. The

System determines whether there is the growth of microorganisms subject to the shining nature change.

2.2.5. Overall dimensions of the Bottles.

Overall dimensions of the Bottles are provided in Table 1.

Dimension	Index
Height, mm	102,0±1,0
Bottom diameter, mm	35,0±0,3
Neck diameter, mm	19,6±0,3
Weight, g, not more than	50±2
Full capacity, ml, not more than	35±1,5
Culture medium capacity, mm, not more than	25±1,5

2.2.6. Culture medium components

Components of culture medium, contained in the Culture medium bottles corresponds to the data of Tables 2 and 3.

Table 2 – Culture medium components

Ingredients list	Culture medium bottles IOHOHA® for aerobe cultivation ■	Antibiotics-neutralized culture medium bottles IOHOHA® for aerobe cultivation for children	Antibiotics-neutralized culture medium bottles IOHOHA® for anaerobe cultivation ■	Antibiotics-neutralized culture medium bottles IOHOHA® for aerobe cultivation ■
Pure water, ml	25,0±1,5	25,0±1,5	25,0±1,5	25,0±1,5
Tryptone, % weight/volume	2,24±0,02	2,22±0,02	2,50±0,03	2,22±0,02
Yeast extract, % weight/volume	0,240±0,005	0,220±0,004	0,250-0,400± (0,005-0,008)	0,220±0,004
Brainheart infusion % weight/volume	0,370±0,007	0,340±0,007	0,100±0,002	0,340±0,007
D(+)-Glucose monohydrate, % weight/volume	0,060±0,001	0,050±0,001	0,060-0,250± (0,001-0,005)	0,050±0,001
Sucrose, % weight/volume	0.080±0,002	0.080±0,002	0,080±0,002	0.080±0,002
D-Fructose, % weight/volume	—	—	0,250±0,005	—
L-Agrinin, % weight/volume	—	—	0,250±0,005	—
Hemin, % weight/volume	—	—	0,00050±0,00003	—
Sodium thioglycolate, % weight/volume	—	—	0,160±0,003	—
Vitamin K1, % weight/volume	—	—	0,000050±0,000003	—

Polyanetholesulfonic acid sodium salt (SPS), % weight/volume	0,060±0,001	0,030±0,001	0,050±0,001	0,050±0,001
Macroporous adsorption resin, g	—	1,6±0,1	1,6±0,1	1,6±0,1

Table 3 – Information on components

Material	Mark	Manufacturer	Country
Tryptone	OXOID	Oxoid Ltd.	United Kingdom
Yeast extract	OXOID	Oxoid Ltd.	United Kingdom
Brainheart infusion	OXOID	Oxoid Ltd.	United Kingdom
D(+)-Glucose monohydrate	Sinopharm	Sinopharm Chemical Reagent Co., Ltd	China
Sucrose	Sinopharm	Sinopharm Chemical Reagent Co., Ltd	China
D-Fructose	AMRESCO	AMRESCO LLC.	USA
L-Agrinin	Sinopharm	Sinopharm Chemical Reagent Co., Ltd	China
Hemin	Sinopharm	Sinopharm Chemical Reagent Co., Ltd	China
Sodium thioglycolate		Shanghai Huixing Biochemical Reagent Co., Ltd.	China
Vitamin K1	Sigma	SIGMA Life Science	Switzerland
Polyanetholesulfonic acid sodium salt (SPS)	Sigma	SIGMA Life Science	Switzerland
Macroporous adsorption resin	DOW	ROHM AND HAAS SHANGHAI CHEMICAL INDUSTRY CO., LTD.	China

2.2.7. Indexes and characteristics requirements

2.2.7.1. Appearance

The Bottle shall be clean in appearance, designations on the label shall be clear. The Bottle shall be leak-proof and shall have no cracks. Leakage and exudation of moisture are not allowed. Culture medium in the Bottle shall be clean, transparent, yellowish with blue resin granules, with no flakes. Polymer sensor shall be blue.

2.2.7.2. Capacity

Culture medium volume in the Bottle amounts to $25 \pm 1,5$ ml.

2.2.7.3. pH value

pH value shall amount to $7,2 \pm 0,2$ at temperature $20-25^{\circ}\text{C}$.

2.2.7.4. Positive results reproducibility

When standard strains are cultured in the Bottles, examination results shall be positive.

2.2.7.4.1. Culture medium sensibility

Culture medium shall ensure standard strains growth not later than 72 hours after the culture at the temperature $35 \pm 1,5^{\circ}\text{C}$. The lowest number of bacteria which shall be determined in the Bottle: 50 ± 10 CFU.

2.2.7.4.2. Growth rate

The minimum time of cultures incubation that shall pass to see a clear and visible to the naked-eye culture growth is 24 hours.

2.2.7.5. Sterility

Examination results shall be negative when standard strains aren't cultured in the Bottles.

2.2.7.6. Consistency

- Once standard strains have been cultured in the Bottles, examination results shall be positive in terms of one month after the period of validity expiration.

- When standard strains aren't cultured in the Bottles examination results shall be negative in terms of one month after the period of validity expiration.

2.2.8. Information on sterility

2.2.8.1. Culture medium sterilization

Ready culture medium shall be sterilized by means of autoclaving in course of 30 minutes at the temperature $(121 \pm 1)^{\circ}\text{C}$, at pressure 0,12MPa. If the culture medium isn't designated for immediate appliance, it shall be stored at temperature $4 \pm 2^{\circ}\text{C}$. After sterilization the culture medium shall be stored for not more than 7 days.

The data is provided in the sterilization process verification report No. XK/TYS-(12)-2015-03.

2.2.8.2. Culture medium Bottles sterilization

Culture medium Bottles exposed to radiation sterilization, which is performed on outsourcing basis by the company Hitrusty Industry and Trade Co., Ltd.

The data is provided in the sterilization process verification report No.HC/FZ-2015-15.

2.3. Composition (delivery completeness) of the medical product

- Culture medium bottles IOHOHa® for aerobe cultivation – 50 pcs. in the box;
- Antibiotic-neutralized culture medium bottles IOHOHa® for aerobe cultivation – 50 pcs. in the box;
- Antibiotic-neutralized culture medium bottles IOHOHa® for anaerobe cultivation – 50 pcs. in the box;
- Antibiotic-neutralized culture medium bottles IOHOHa® for aerobe cultivation for children – 50 pcs. in the box;
- Instruction for use;
- Certificate.

2.4. Warnings, safety measures and restrictions during medical product use

2.4.1. Warnings and safety measures

- Learn instruction for use before working with the Bottles;
- The Bottle with evidence of damage or leakage aren't allowed to be used;
- The culture medium bottles are designated for single use;
- Check whether the text, designations and marking on the label are clear before using the Bottles.

2.4.2. Restrictions for use

- Bacteria with special nutritious demands, defined by slow rate of growth, in a sample can result in its insufficient growth in the Bottle. There's a possibility of false negative result.
- Results mostly depend on blood samples collection conditions, so one negative result doesn't accurately determine infection absence. It's required to comply with the cultivation procedure.

2.5. Package, marking and warning designations description

2.5.1. Package

Culture medium is packed in the Bottle which ensures protection of the medium from mechanical and climatic factors during transportation and package.

Culture medium bottles are put in transportation package, which is a carton box with dimensions: 393x195x115 mm, and high strength corrugate paper is put inside in order to prevent the Bottles consistency damage. Carton box design allows to exclude culture medium bottles movability during transportation and to protect products from damage.

2.5.2. Marking


2.5.2.1. Bottle marking

Culture medium bottle label contains the following information:

- manufacturer's logo,
- culture medium volume,
- culture medium name and designation,
- sample volume,
- sign "Medical product for in vitro diagnostics" (corresponding symbol is applied);
- storage conditions:
 - sign "Protect from moisture" (corresponding symbol is applied);
 - sign "Temperature range" (corresponding symbol is applied);
 - sign "Prevent exposure to sunlight" (corresponding symbol is applied);
- designation;
- product standard;
- license number;
- number of registration license in the Russian Federation;
- manufacturer's name;
- manufacturer's country;
- bar code;
- detachable bar code;
- batch number;
- production date;

- validity period with the sign “expire date”;
- place for information on patient:
 - last name
 - time
 - number

Indication of applied symbols:

	Medical product for <i>in vitro</i> diagnostics
	Sterile
	Protect from moisture
	Prevent exposure to sunlight
	Temperature range

2.5.2.2. Identification by colour

Culture medium bottles differ from each other by label and removable cap colour (images 5a, 5b, 5c, 5d, 5e) depending on its designation:

- Culture medium bottle **ЮHOHA®** for aerobe cultivation – red:

Image 5a. Bottle label sample

- Antibiotic-neutralized culture medium bottle ЮНОНА® for aerobe cultivation – blue:

Image 5b. Bottle label sample

- Antibiotic-neutralized culture medium bottle ЮНОНА® for anaerobe cultivation – cadmium orange:


Image 5c. Bottle label sample

- Antibiotic-neutralized culture medium bottle ЮНОНА® for aerobe cultivation for children – pink:


Image 5d. Bottle label sample

2.5.2.3. Identification by bar code

Culture medium bottles differ from each other by identification code, contained in bar code, depending on its designation:

Name	Sample volume	Culture medium volume	Designation	Bar code/ Identification number	Photo
Culture medium bottle ЮНОНА® for aerobe	5-10 ml	25 ml	biological samples of	XBZ	

cultivation			patients who didn't take antibiotics		
Antibiotic-neutralized culture medium bottle ЮHOHA® for aerobe cultivation	5-10 ml	25 ml	biological samples of patients who took antibiotics	XCR	
Antibiotic-neutralized culture medium bottle ЮHOHA® for anaerobe cultivation	2-10 ml	25 ml	biological samples of patients who took antibiotics	XYY	
Antibiotic-neutralized culture medium bottle ЮHOHA® for aerobe cultivation for children	2-10 ml	25 ml	biological samples of children and other samples difficult to collect	XET	

2.5.2.4. Transport package marking:

- manufacturer logo;
- manufacturer name;
- handling signs:
 - top (corresponding symbol is applied),
 - protect from moisture (corresponding symbol is applied),
 - regular pile legs number limit (corresponding symbol is applied),
 - do not roll (corresponding symbol is applied),
 - designation “Medical device for in vitro diagnostics” (corresponding symbol is applied);
 - information on sterility (corresponding symbol is applied);
 - culture medium name and designation;
 - batch number;
 - manufacture date;
 - validity period;
 - storage temperature;
 - number of pieces in a box;
 - manufacture address;

- product standard;
- license number;
- number of registration license in the Russian Federation ,
- manufacturer address and telephone;
- authorized representative name, address and telephone.
- market circulation mark;
- weight netto/brutto.

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


Image 5e. Ssample transport package marking

Applied symbols designations:

	Top
---	-----

	Protect from moisture
	Regular pile legs number limit
	Do not roll
	Medical product for <i>in vitro</i> diagnostics
	Sterile
	Market circulation mark

3. Medical product operation

3.1. Adjustment, assemblage, setting up, gauging and other activities, required to put the product into operation.

3.1.1. Adjustment and assemblage aren't required for using the Bottles as intended.

3.1.2. The Bottles are designated for joint operation with Automated bacteriological System **ЮHOHA®** LABSTAR. See the System operation conditions in the operation manual for the medical product "Automated bacteriological System **ЮHOHA®** LABSTAR, models: **ЮHOHA®** LABSTAR 50, **ЮHOHA®** LABSTAR 100, with accessorises».

3.2. Operation procedure

3.2.1. General regulations for collecting biological samples for the examination.

- It's reasonable to collect biological material before starting antibacterial therapy.
- Material for bacteriologic examination shall be collected directly from the infection pocket or clinically significant biological material shall be examined.
- It's required to respect antiseptics and avoid biological material contamination with foreign bacterial population.
- Material amount shall be sufficient for examination realization.
- It is not allowed to use the Bottles with evidence of damages and leaks.
- Culture medium bottles are designated for the single use.
- Make sure that the text, signs and marking on the Bottles labels are clear before using the Bottles.
- If anaerobic infection is suspected, a sample should be collected in two vials simultaneously - with culture media for aerobic and anaerobic cultivation.
- Do not use expired medium
- Culture medium transparency shall be visually determined before using the Bottles, and any turbidity shall evidence its invalidity.
- Before sowing the biological material in the bottle of the medium should be warmed to a temperature of 18-30°C.

3.2.2. Blood examination

- If temperature increases the blood shall be collected from the peripheral vein, antiseptics regulations and individual protection measures, required while working with blood, shall be respected (rubber gloves shall be used).
- Sterile syringes and systems for single use are required for the blood collection.
- Exhaustion is created in the Bottle, and it simplifies patients' blood collection.

- Blood volume from 5 ml to 10 ml for adults for aerobic bottles, from 2 ml to 10 ml for anaerobic bottles and from 2 ml to 10 ml for children (the label on the label of the vial allows you to estimate the volume of blood taken). For newborns, the amount of blood taken should not exceed 1% of the total blood volume. Optimal pathogen isolation is achieved by adding the maximum sample volume.

- Blood shall be cultures in culture mediums straight after its collection at the patient's side or in the medical treatment room.

- If the patient is suspected of bacterial endocarditis, three microbiological tests are performed to isolate the pathogens, and blood is drawn at least once per hour for 24 hours.

- If the patient has been treated with antibiotics for one to two weeks prior to the blood culture test, then two to three microbiological blood tests should be performed within one day.

- Blood shall not be collected with intravascular catheters excepts if catheter origin infection is suspected.

3.2.2.1. Blood collection procedure and culture

- Culture medium bottles, sterilizing materials and other materials shall be prepared (image 6, drawing 1).

- Culture medium bottle shall have fluid volume scale, where 5-10 ml blood filling volume shall be designated (image 6, drawing 2).

- Plastic protection cap shall be removed from the Bottle, the Culture medium bottle shall be sterilized with the cotton wool stick/wipe, soaked in 70% ethyl hydroxide, avoid using iodine. After sterilizing alcohol wipe/stick shall be applied to the Bottle for sterility maintenance (image 6, drawing 3).

- Skin disinfection: the skin shall be disinfected with the stick/wipe, soaked in 70% ethyl hydroxide and left to dry for 30 seconds. Then a center point for the puncture shall be marked a with iodine which is required for disinfection. Then a cycle of 2.5-3 cm radius shall be disinfected (1-2% iodine tincture for 30 seconds) (image 6, drawing 5).

- Iodine shall be removed with the wipe/cotton wool stick soaked in 70-75% ethanol solution (image 6, drawing 6).

- It's recommended to use a pair of aerobic and anaerobic bottle for adult patients

- It's recommended to use the acus for vacuum blood collection ((butterfly needles with Luer adapter, sterile, 21-23G) for collecting blood directly to the Bottle. It excludes personnel contact with the blood, protects personnel from spontaneous acus wound, which can lead to hematogenous infection (image 6, drawing 7).

- When using a blood sampling system, culture is done first in an aerobic vial and then in an anaerobic vial. (Do not remove the needle from the patient's vein during blood collection. After the end of the first extract needle from vial, then from veins patient). When using a syringe, first inject blood into an anaerobic vial and then into an aerobic one (avoid air entering the anaerobic bottle.)

- The Bottle shall be turned up-down carefully (image 6, drawing 8).


Image 6 Patient's blood collection procedure

- After the procedure, iodine is removed from the skin with a cloth with 70% ethyl alcohol and the puncture area is closed with a plaster.

- The Bottles shall be marked with the information on the patient and material collection date. They shall be stored at indoor temperature (18-30°C) in a place protected from the light (not in the fridge) before their transporting to the laboratory.

- It's required to remove half of a label, which contains bar code, from the Bottle (it is impossible to remove the second part) and stick it to the laboratory blood examination assignment module.

When removing the label it is required to check whether the information on the assignment and the Bottle coincides!

⚠ It's required to send samples to the laboratory for examination as soon as possible. It will allow to increase the possibility of pathogens detection. It is required to start sample examination not later than in 24 hours after collection!

⚠ If samples processing cannot be started, they shall be stored at indoor temperature, do not put them in the fridge! In winter time culture bottles should be delivered to the laboratory in a thermal container with a temperature of 18-35°C.

3.2.3. Spinal fluid examination

- Cerebrospinal fluid (CSF) is obtained with lumbar puncture, rarely with puncture of the lateral ventricles of the brain, preferably before taking antibiotics

- The puncture area is treated with 70% ethyl alcohol, then 1-2% iodine tincture for 30 seconds. The first portions of liquor have the highest potential for contamination by microorganisms from the skin and should be discarded.

- Immediately after puncture, not less than 5 ml and not more than 10 ml of the resulting liquor should be aseptically introduced into an aerobic vial (in certain clinical situations, an anaerobic vial is required secondarily).

If the volume of liquor is small from 2 to 5 ml, as well as liquor from children - you should use aerobic bottles for children.

- Large volumes (5-10 ml) are introduced into bottles with an aerobic medium.

- Recommended volume for anaerobic culture medium bottles is 2-10 ml.

- Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of pathogens.

- During storage and transportation, excessive overheating and especially cooling should be avoided (a thermos or thermal container with a temperature of 30-35°C can be used). Cooling the liquor below 30°C leads to the death of meningococci.

- Storage and transportation time should not be more than 8 hours.

- After sowing liquor in the bottle, you must tear off half of the label with a bar code from the bottle (the second half of the label can not be torn off) and glue it to the form of the direction for laboratory testing of liquor.

 When tearing off the label with a bar code, it is necessary to ensure that the information on the direction and on the bottle matches!

3.2.4. Examination of body fluids from normally sterile cavities:

- Peritoneal, synovial, pleural, pericardiac fluids can be examined.
- Fluids collected by means of the puncture and aspiration shall be used for bacteriological examination.
 - The puncture area is treated with 70% ethyl alcohol, then 1-2% iodine tincture for 30 seconds
 - For adults, it is recommended to take for each aerobic vial from 5 ml to 10 ml, for each anaerobic vial from 2 ml to 10 ml of the test liquid, for children, as well as in the case of restrictions on the volume of the collected sample for each vial you need from 2 ml to 10 ml of the test liquid
 - The rate of growth of pathogens is directly proportional to the volume of biomaterial. Optimal pathogen detection is achieved with maximum sample volume
 - Immediately after the puncture, the biological fluid should be aseptically introduced into the aerobic vial (in certain clinical situations, an anaerobic vial is required secondarily).
 - Air bubbles, entered in a syringe; shall be removed.
 - Minimum fluid volume for bacteria discharge is 2-5 ml, minimum fluid volume for fungus discharge is 10 ml.
 - Excessive fluid and purulence shall be transported in sterile containers with screw cap.
 - Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of microorganisms
 - Store at room temperature (18-30°C) for a maximum of 24 hours.
 - After sowing sample in the bottle, you must tear off half of the label with a bar code from the bottle (the second half of the label can not be torn off) and glue it to the form of the direction for laboratory testing of liquor.

 When tearing off the label with a bar code, it is necessary to ensure that the information on the direction and on the bottle matches!

3.2.5. Cultivation

3.2.5.1. The Bottles with cultured biological material shall be put in the main block of **ЮЮЮЮЮ**[®] LABSTAR automated bacteriological system.

3.2.5.2. If System alarm is tripped, software interface shall be checked:

- If the position on the screen flashes red, it means positive result. There is growth of microorganisms in the Bottle. It is required to remove the Bottle from the System, execute reculture and send the sample for additional examination. More than 90% of positive samples are detected within 24 hours.

- If the position on the screen turns green when the maximum cultivation time is set to 120 h, it means a negative result.. There is no growth of microorganisms in the Bottle or there are inviable microorganisms.

Additional identification and drug sensitivity can be executed for colonies taken from pure growth.

3.2.5.3. Received results interpretation

- Positive results:
 - There is growth of microorganisms in the blood.
 - Sample contamination during collection can lead to a false positive result.
 - Execute culture for confirmation.
- Negative results:
 - There is no growth of microorganisms in the blood.
 - There are bacteria with special nutritious demands, defined by slow rate of growth.
 - There is a small amount of bacteria in a sample.
 - Antibiotics intake can cause false negative result.
 - Storage condition violation can influence efficiency of microorganisms detection and even lead to a false negative result.

3.3. Quality control

The certificate of quality of the cultural media bottles is provided by the manufacturer with each batch of the product.

In-laboratory quality control of the media bottles is not required. However, if desired, laboratories can perform media quality control using standard strains (see Table 4) and the LABSTAR automated bacteriological system.

Table 4. Reference ATCC strains for examination

Bottles with aerobic cultures	Bottles with anaerobic cultures
<i>Neisseria meningitidis</i> (ATCC13090)	<i>Clostridium histolyticum</i> (ATCC19401)
<i>Haemophilus influenzae</i> (ATCC19418)	<i>Clostridium perfringens</i> (ATCC13124)
<i>Streptococcus pneumoniae</i> (ATCC6305)	<i>Bacteroides fragilis</i> (ATCC25285)
<i>Streptococcus pyogenes</i> (ATCC19615)	<i>Bacteroides vulgatus</i> (ATCC8482)
<i>Pseudomonas aeruginosa</i> (ATCC27853)	<i>Streptococcus pneumoniae</i> (ATCC6305)
<i>Candida albicans</i> (ATCC 18804)	<i>Escherichia coli</i> (ATCC25922)
<i>Staphylococcus aureus</i> (ATCC25923)	<i>Staphylococcus aureus</i> (ATCC25923)
<i>Escherichia coli</i> (ATCC25922)	
<i>Alcaligenes faecalis</i> (ATCC8750)	

Positive and negative control must be exercised simultaneously. For positive control, prepare an inoculate with microorganisms of one or more of the types specified in Table 4. as follows:

- Using a sterile physiological 0.9% solution to prepare inoculate 0.5 MCF (0.5 units of turbidity according to the Mcfarland standard) of the reference strains of table 4 (tube number 0-0.5 u = $1.5 \cdot 10^8$ CFU / ml);
- Take another 4 tubes of 10 ml sterile physiological 0.9% solution-respectively tubes number 1,2,3,4);
- Take 100 μ l from test tube No. 0 and add to test tube No. 1. Mix thoroughly (the concentration of bacteria in test tube No. 1 will be $1.5 \cdot 10^6$ CFU/ml) ;
- Take 100 μ l from test tube No.1 and add to test tube No.2. Mix thoroughly (the concentration of bacteria in test tube No. 2 will be $1.5 \cdot 10^4$ CFU/ml) ;
- Take 100 μ l from test tube No. 2 and add to test tube No. 3. Mix thoroughly (the concentration of bacteria in test tube No. 3 will be $1.5 \cdot 10^2$ CFU/ml);
- Take 5 ml from test tube No. 3 and add to No. 4. Mix thoroughly (the concentration of bacteria in test tube No. 4 will be 50 CFU/ml).

Syringe 1 ml of bacterial suspension 50 CFU/ml into bottle of medium and place them in the system.

Growth of microorganisms in the seeded medium is usually observed within 72 hours, more than 90% of positive samples are detected within 24 hours. In this case, the warning system is activated, and the position on the system screen with a positive bottle lights up in red.

In the absence of seeding, the result should be negative. The system is activated after 5 days (120 h is the default cultivation time). The medium in the bottle should be clean and transparent.

3.4. Medical product application risks, counterindications, expected and predictable side effects, connected to medical product use as intended.

3.4.1. The Bottles are not accepted for examination, if:

- Information of the Bottle does not coincide with the analysis conduction assignment.

- The Bottles have cracks or visible impurity.
- 24 hours passed after collection.
- The Bottle type does not correspond to examination type.
- The Bottle validity period expired.
- Material volume in the Bottle is not sufficient for examination.

3.4.2. If there are bacteria with special nutritious demands, defined by slow rate of growth, it can lead to their insufficient growth in the Bottle. There is a possibility of a false negative result.

3.4.3. The results mostly depend on collection conditions, so one negative result does not determine absence of infection. It is required to respect cultivation procedure strongly.

4. Information on maintenance and current repairs

Maintenance is not foreseen. The Bottle is a single use product.

5. Storage and transportation

5.1. Transportation

- The Bottles transportation shall be performed at the temperature 4-30°C. It is required to take corresponding measures for necessary temperature maintenance in conditions of cold winter or hot summer.

- It is required to protect the Bottles from low temperatures and also from high temperature continuance.

- It is required to exercise caution, avoid pressure, keep package clean, not to open the package, protect it from water, atmospherical condensation and sunlight impact.

5.2. Storage

- It is required to store packages with the Bottles in a dry place at the temperature 4-30°C.

- It is required to protect it from direct sunlight.
- The Bottles period of validity is 12 months.

6. Warranty liabilities

- The Bottles period of validity is 12 months from their production date.
- The manufacturer guarantees the Bottles compliance with applicable standard and manufacturer documentation requirements if conditions of transportation, storage and operation are met.
- The manufacturer is not liable for any damage or defect done during the Bottles use until his guilt has been proved.
- The manufacturer is not liable for any damage or defect, caused but repeated use of the Bottles.
- The manufacturer is not liable for any damage or defect, caused by the use of the Bottles after expiration of the validity period, specified on the package.

7. Medical product disposal and liquidation procedure.

The product must be disposed of in accordance with the current legislation of the Russian Federation at the time of disposal.

Samples may contain pathogens, including hepatitis B and human immunodeficiency virus (HIV), therefore, the media belong to class B - "Epidemiologically hazardous waste" and require appropriate precautions for their disposal according to SanPiN 2.1.7.2790-10.

Used media shall be disinfected. If the laboratory is equipped with an autoclave, the disinfection is performed by autoclaving (exposure to water vapor under pressure) at a temperature of 121°C and a pressure of 1.2 ATM. 20 minutes. Vials with the remnants of blood or other biological samples, closed lids, harvest in sustainable for leaks fluid capacity or plastic packages, suitable for physical method disinfection and housed in an autoclave. After autoclaving, the waste is transferred to yellow bags labeled "class B Medical waste". The person responsible for the collection of waste in this medical unit ties the bag or closes it using tie-tags or other devices and indicates the date, time, Name of the person responsible on the label. Packages are removed to the temporary storage of waste for subsequent removal and disposal. If the laboratory is not equipped with an autoclave, disinfection is carried out using other methods established by local rules and regulations.

Medium that have become unusable because of a violation of the integrity of the packaging (bottles) as well as expired can be disposed of in accordance with SanPiN 2.1.7.2790-10 as medical waste class B "Epidemiologically hazardous waste".

Cardboard boxes and liners are disposed of in accordance with SanPiN 2.1.7.2790-10 as class A waste "Epidemiologically safe waste, in composition close to solid household waste."

Appendix A

The information about the examination with the use of American typical culture collection (ATCC) strains.

Strains	Strain code	Conditions		Pathogenic microorganism detection method	Results
<i>Neisseria meningitidis</i>	ATCC13090	Temperature (°C)	Time (h)	The culture detected by System	Triggering the alerting system
<i>Haemophilus influenzae</i>	ATCC19418	35	<72	The culture detected by System	Triggering the alerting system
<i>Streptococcus pneumoniae</i>	ATCC6305	35	<72	The culture detected by System	Triggering the alerting system
<i>Streptococcus pyogenes</i>	ATCC19615	35	<72	The culture detected by System	Triggering the alerting system
<i>Pseudomonas aeruginosa</i>	ATCC27853	35	<72	The culture detected by System	Triggering the alerting system
<i>Candida albicans</i>	ATCC18804	35	<72	The culture detected by System	Triggering the alerting system
<i>Staphylococcus aureus</i>	ATCC25923	35	<72	The culture detected by System	Triggering the alerting system
<i>Escherichia coli</i>	ATCC25922	35	<72	The culture detected by System	Triggering the alerting system
<i>Alcaligenes faecalis</i>	ATCC8750	35	<72	The culture detected by System	Triggering the alerting system
<i>Clostridium histolyticum</i>	ATCC19401	35	<72	The culture detected by System	Triggering the alerting system
<i>Clostridium perfringens</i>	ATCC13124	35	<72	The culture detected by System	Triggering the alerting system
<i>Bacteroides fragilis</i>	ATCC25285	35	<72	The culture detected by System	Triggering the alerting system
<i>Bacteroides vulgatus</i>	ATCC8482	35	<72	The culture detected by System	Triggering the alerting system

Appendix B

Diseases, caused by microorganisms

Aerobic bacteria:

➤ **Neisseria meningitidis** are pathogens which cause children's and adults' cerebral fever, bacteraemia and other serious infections.

➤ **Haemophilus influenzae** can cause acute suppurative inflammation, bacteraemia, sepsis. Existing capsular microbes, especially b-type capsular microbes may cause primary infection. Non-capsular bacteria are also pathogens and they mostly cause secondary infection. Their virulence is mostly caused by bacterial endotoxins, which can cause such patients' secondary diseases as flu, measles, whooping cough, bronchitis. After infection anticapsular specific antibodies appear, and resistant immunity forms.

➤ **Streptococcus pneumoniae** may cause sudden strong fever, depressed state of mind, impatience, anxiety, headache, cough, laborious breathing, otalgia, dysacusis, purulent outflow and other symptoms.

➤ **Streptococcus pyogenes** is mainly the cause of acute pharyngitis, acute tonsillitis, and it also can cause pulmonary infection, scarlet fever, skin covering and soft tissues infection, and it also can cause generalized infection. These bacteria are also indirect cause of the diseases of allergic Bouillaud's disease and acute glomerulonephritis.

➤ **Pseudomonas aeruginosa** is a relative pathogen, and it is the natural barrier and does not cause pathogen changes if skin covering is not damaged, but if skin covering is damaged, it can cause the following diseases: septic fever, respiratory infections, endocarditis, urinary tract infections, brain infections and so on.

➤ **Candida albicans** can be detected in different parts of the human body and it can cause the following diseases: cutaneous candida infection (mainly it appears in skin folds, erythema, intertrigo, skin, sometimes the skin covers with whitish plaque and cracks, affected areas with vesicles); mucous membrane candida infection (in most cases it causes stomatitis, angular cheilitis, vaginitis, mucous membrane surface covers with white curdled plaque with large grains. After curdled-white formations removal, superficial painful erosion is left); internal organs and cerebrospinal nervous system candidiasis (includes pneumonia, gastroenteritis, cerebral fever, encephalitis and so on, and in rare cases it includes septic fever).

➤ **Staphylococcus aureus** may cause allocated purulence infection, and also it can cause pneumonia, pseudomembranous colitis, pericarditis and even septical fever, pyemia and other generalized infection.

➤ **Escherichia coli** diseases, caused by escherichia coli mainly cause modified bronchoalveolitis in lower lobes of both lungs. On the 6th day of disease development patients often have signs of pulmonary abscess, pleural fluid and even pyothorax. There is lower possibility of inflammation in mucous membrane of trachea and bronchial tubes. It may be so because mostly when a person suffers from pneumonia caused by escherichia coli, the contamination happen in a haematogenous way. There is a common amount of serous liquid in bronchic cells and a common amount of mononuclear cells. On early stages of disease development a large amount of erythrocytes appears in exudate, on later stages of disease neutrophils and macrophages are detected. Cells walls enlargement and visible necrotic damages can be observed. In some cases escherichia coli can cause cholecystitis, pyelonephritis, and cerebral fever, and other diseases.

➤ **Alcaligenes faecalis** causes enteritis, urinary tract infections.

Anaerobic bacteria:

➤ **Clostridium histolyticum**

➤ **Clostridium perfringens**

➤ **Bacteroides fragilis**

➤ **Bacteroides vulgatus**

Anaerobic infection is a hard toxic wound fever, caused by anaerobic microorganisms. Endogenous infection develops as a result of anaerobic migration inside the body. Infections, caused by anaerobic bacteria, are usually provoked by mixed flora. If a negative result of bacterial culture cultivation is detected, the possibility of anaerobic infection presence is not excluded. If a positive result of bacterial culture cultivation is detected, the possibility of mixed anaerobic infection is also not excluded.

Clinical criteria of anaerobic infections presence include the following:

➤ Active gas generation, tissues edema and necrosis, subcutaneous crackling rale are present in infected tissues. These are symptoms, which are caused by **Clostridium perfringens**.

➤ The infection affects mucous membrane, affected area periphery, oral cavity, intestinal canal, pharyngonasal cavity, vagina and other mucous membranes easily, here large number of anaerobic parasites are concentrated. If these parts of

the body and their periphery are damaged, they are highly exposed to anaerobic infection;

- Such deep traumas, as gunshot wounds and animal bites create supportive environment for anaerobic infections;

- If the secretion has stinking odor or dark flesh red colour, and it has red fluorescence at UV light, it corresponds to anaerobic infection. If the secretion or the pus have "grey granules", it corresponds to antimycotic infection;

- If microbes are detected in the swipe sample of patient's secretion microscopically after Gram coloration, but after cultivation or bacteria cultivation in fluids or in deep layers of semisolid culture medium a negative result is obtained, anaerobic infection can be detected;

- In case of long-term aminoglycoside antibiotics intake and their inefficiency, the anaerobic infection is possible;

- Infection can be caused by abortions and gastrointestinal surgical operations;

- It shall be considered that if the bacterial reaction is negative, and bacterial endocarditis, detected in the blood culture, is complicated with septic thrombophlebitis and accompanied by yellow plague and so on, it can state anaerobic infection.

Thus, the use of culture media bottles for microbiological diagnosis of infectious diseases should be combined with a number of additional bacteriological and clinical studies and taking into account the history of patients.

Appendix C

List of applicable national standards to ensure the safety, efficacy and quality of medical products

Medical product "Culture medium Bottles ЮHOHA® for Automated bacteriological Systems ЮHOHA® LABSTAR in models" complies with international and national standards of the People's Republic of China for products:

- GB/T2828.1-2012 Sampling procedure for quality control-Part 1: Index of sampling schemes according to acceptable quality limits (AQL) for sequential batch inspection;
- GB/T2829-2002 Procedure for sampling and the frequency of periodic quality control (Applicable for the stability control of the process);
- GB/T19973.2-2005 Sterilization of medical products-Microbiological methods-Part 2: Sterility Tests carried out during the validation of the sterilization process;
- ISO 11737-1: 2006 Sterilization of medical products. Microbiological methods - Part 1: Determination of a population of microorganisms on products;
- ISO 11137-1: 2006 Sterilization of medical products-Radiation sterilization-Part 1: Requirements for the development, validation and ongoing monitoring of the sterilization process of medical products;
- ISO 11137-2: 2012 Sterilization of medical products-Radiation sterilization-Part 2: Establishment of sterilizing dose.
- YY0033-2000 Good manufacturing practices for sterile medical products;
- YY0466.1-2009 Medical Products-Symbols on labels of medical products and the provided information, Part 1: General requirements
- YY / T 0656-2008 Automated Blood Culture System;
- Principles and procedure for handling blood cultures. Approved methodical instruction;
- National regulations on the organization and operation of clinical laboratories (Third edition);
- Chinese Pharmacopoeia.

Medical product "Culture medium Bottles **ЮHOHA®** for Automated bacteriological Systems **ЮHOHA®** LABSTAR in models" complies with international and national standards of the Russian Federation for products:

- GOST R 51352-2013 «In vitro diagnostic medical devices. Test methods»;
- GOST R 51088-2013 «In vitro diagnostic medical devices. Reagents, kits, the test - systems, control materials, culture media. Requirements to devices, and to supporting documentation»;
- GOST R ISO 23640-2015 «In vitro medical devices. Evaluation of stability of in vitro diagnostic reagents»;
- GOST R EN 12322-2010 « In vitro diagnostic medical devices. Culture media for microbiology. Performance criteria for culture media»;
- GOST R ISO 18113-1-2015 «In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 1. Terms, definitions and general requirements»;
- GOST R ISO 18113-1-2015 «In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 2. In vitro diagnostic reagents for professional use»;
- GOST R ISO 15223-1-2014 «Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements»;
- GOST R 56894-2016 «Summary technical documentation for demonstrating conformity to the essential principles of safety and performance of in vitro diagnostic medical devices»;
- SanPiN 2.1.7.2790-10 «Sanitary and epidemiological requirements for the treatment of medical waste».

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