Operational Documentation For Medical Device

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Name and address of Production site:

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Approved by:

Full name: Stéphanie HOMOLLE
Position: International Regulatory Affairs Specialist
Company Name: bioMérieux S.A.
Signature: *[Handwritten Signature]*
Date: *January, 20th, 2014*
Stamp:

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000 ЦРММ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Content of Technical Documentation according to Appendix G of the “Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices”		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Name of the device (IVD MD)	X	Yes	VITEK MS-CHCA Top of the first page of Package Insert (PI)
	Purpose of the IVD MD	X	Yes	Top of the first page of PI VITEK MS-CHCA matrix is an auxiliary reagent for <i>in vitro</i> diagnostic and microbiological control. It is intended to be used for the preparation of samples for further analysis by means of MALDI-TOF (Matrix-Assisted Laser Desorption/Ionization – Time Of Flight), exclusively on VITEK MS mass spectrometer, manufactured by bioMérieux S.A., France (registration certificate No. FSZ 2012/12091 from 05.05.2012). A set of reagents (including VITEK MS-CHCA matrix) is used to conduct a compulsory stage of preparation of microorganisms sample. Further in the analysis process, the matrix absorbs the laser light and vaporizes, along with the sample, in the process gaining an electrical charge (ionization).
	Indication that the IVD MD is intended for use for clinical laboratory diagnostics	X	Yes	Section: Warnings and Precautions of PI
	Designated user	X	Yes	Section: Warnings and Precautions of PI

Content of Technical Documentation according to Appendix G of the “Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices”		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Test principle		N/A	Not applicable in this case: VITEK MS-CHCA matrix is an auxiliary reagent for <i>in vitro</i> diagnostic and microbiological control. It is intended to be used for the preparation of samples for further analysis by means of MALDI-TOF (Matrix-Assisted Laser Desorption/Ionization – Time Of Flight), exclusively on VITEK MS mass spectrometer, manufactured by bioMérieux S.A., France (registration certificate No. FSZ 2012/12091 from 05.05.2012). A set of reagents (including VITEK MS-CHCA matrix) is used to conduct a compulsory stage of preparation of microorganisms sample. Further in the analysis process, the matrix absorbs the laser light and vaporizes, along with the sample, in the process gaining an electrical charge (ionization). For more information, see: • Section 4 of User Manual - 4501-2233 - C - en - NA - VITEK MS Clinical Use – Workflow • Section 4 of User Manual - 4501-2233 - C - ru - NA - VITEK MS Clinical Use – Workflow
	Description of reagents, calibrators and control materials.	X	Yes	Section: Content of the kit of PI The kit does not contain any calibrators and control materials.

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	- List of the delivered materials making part of the IVD MD - List of accessories - List of special materials required to conduct an analysis, but not delivered with the MD set		N/A	Not applicable in this case: The kit does not contain any accessories nor need special materials
	Indication of any special storage conditions (e.g. temperature, light, humidity, etc.) and/or special conditions of handling the IVD MD	X	Yes	Section: Storage conditions of PI
	Stability parameters (storage conditions, shelf life after the first opening, storage conditions and stability of working solution, where applicable)	X	Yes	Section: Storage conditions of PI
	For IVD MD delivered as a sterile product - indication to its sterile condition, sterilization method and instructions in case of damaging sterile packaging.		N/A	Non-sterile
	Information regarding warnings, precautionary measures and restrictions on the use of the IVD MD.	X	Yes	Section: Warnings and Precautions of PI

25020 411071
MDA 100

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	- List of the delivered materials making part of the IVD MD - List of accessories - List of special materials required to conduct an analysis, but not delivered with the MD set		N/A	Not applicable in this case: The kit does not contain any accessories nor need special materials
	Indication of any special storage conditions (e.g. temperature, light, humidity, etc.) and/or special conditions of handling the IVD MD	X	Yes	Section: Storage conditions of PI
	Stability parameters (storage conditions, shelf life after the first opening, storage conditions and stability of working solution, where applicable)	X	Yes	Section: Storage conditions of PI
	For IVD MD delivered as a sterile product - indication to its sterile condition, sterilization method and instructions in case of damaging sterile packaging.		N/A	Non-sterile
	Information regarding warnings, precautionary measures and restrictions on the use of the IVD MD.	X	Yes	Section: Warnings and Precautions of PI

29230 411071
NOV 100

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
	Warnings and/or precautionary measures associated with potentially infectious materials contained in the IVD MD.	X	Yes	- Section: Warnings and Precautions of PI <i>Note: VITEK MS-CHCA does not contain any potentially infectious materials. However all specimens, microbial cultures and inoculated products to be analyzed on VITEK MS instrument should be considered infectious and handled appropriately.</i>
	As applicable: - special requirements for the premises - special training - special qualification of the expected user of the IVD MD.	X	Yes	Section: Warnings and Precautions of PI
	- Conditions for collection, processing and preparation of samples - Information about stability of the analyzed samples		N/A	Not applicable in this case: the VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization.
	Information required to verify whether the IVD MD is installed correctly and is ready to for safe operation	X	Yes	Section: Warnings and Precautions of PI

29.03.2014 11:10:12
NDN 100

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Recommendations regarding the quality control procedures, where applicable.	X	Yes	Section: Quality Control of PI
	Information about metrological traceability of the values assigned to the calibrators and the materials for control of correctness, including an indication to the applicable reference materials and/or reference measuring procedures of higher levels.		N/A	Not applicable in this case: The kit does not contain any calibrators and control materials.
	Analysis procedure, including calculations and interpretations of the results, and, where applicable, information about the feasibility of confirmation tests.		N/A	Not applicable in this case: • VITEK MS-CHCA does not have any measuring function. • VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization.
	Analytical parameters, such as sensitivity, specificity, correctness, repeatability, reproducibility, detection limits and measurement range, including information regarding control of known interfering agents, restrictions of the method and use of available reference materials and measurement methods.		N/A	Not applicable in this case: the VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization.

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Clinical efficacy parameters, such as diagnostic sensitivity and diagnostic specificity, as applicable.		N/A	Not applicable in this case: the VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization.
	Biological reference interval, where applicable.		N/A	Not applicable in this case: the VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization.
	Information about interfering substances or restrictions associated with the sample, which can impact the study result.		N/A	Not applicable in this case: the VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization.
	Warnings and/or precautionary measures associated with the safe disposal of the IVD MD and accessories, if any	X	Yes	Section: Waste disposal of PI

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Name of the manufacturer and/or its authorized representative, as well as their addresses and contact details, from which technical advice and support can be obtained.	X	Yes	Bottom of last page of PI
	Instruction of the user of the necessity to furnish a notice containing the respective information to the Federal Service on Surveillance in Healthcare of the Russian Federation according to the applicable laws in case of detection of any side effects which are not specified in the instruction for use or the operation manual of the medical device, the adverse reactions during its use, the peculiarities of interaction of medical devices one with another, the facts and circumstances posing threat to the life and health of citizens and healthcare professional in the process of application and operation of the medical devices.	X	No	See bottom of last page of Technical and Operational Documentation for VITEK MS-CHCA.
Section 2: Additional instructions of use for the IVD MD intended for self-testing or testing near the sickbed	Detailed information about the analysis procedure, including preparation of reagents, taking/preparation of a sample, and information regarding the procedure for conducting the test and calculating the results.		N/A	Not applicable in this case: VITEK MS-CHCA matrix is not intended for self-testing nor for near patient testing.

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 2: Additional instructions of use for the IVD MD intended for self-testing or testing near the sickbed	- Recommendations regarding the user's actions in case of a positive, negative or indefinite result; - Information about the test restrictions and the possibility of obtaining falsely positive or falsely negative results, as well as information regarding any factors which may influence the test result (e.g. age, sex, menstrual periods, infections, exercises, diets or medicine administration).		N/A	Not applicable in this case: VITEK MS-CHCA matrix is not intended for self-testing nor for near patient testing.
	For IVD MD intended for self-testing - an indication to the fact that the user should not make medical decisions without the prior consultation with the healthcare professional.		N/A	Not applicable in this case: VITEK MS-CHCA matrix is not intended for self-testing nor for near patient testing.

On issues related to the quality of medical products "VITEK MS-CHCA", and in case of detecting any side effects specified neither in the instruction for use nor the operating instructions on the medical devices, of adverse reactions when used, interaction of medical devices, as well as facts and circumstances endangering the life and health of citizens and medical professionals during the use and operation of medical devices please contact Authorized representative of the manufacturer: **LLC "bioMerieux RUSS"**

Address: Russia
115114, Moscow
Derbenevskaya Street
Building 20
Structure 11

Tel./fax: (495)221 10 79
e-mail: office.ru@eu.biomerieux.com
Website: www.biomerieux-russia.com

The reports on possible side effects not specified in the instruction for use or the operating instructions on the medical devices, of adverse reactions when used, interaction of medical devices, as well as facts and circumstances endangering the life and health of citizens and medical professionals during the use and operation of medical devices will be sent to the Federal Service for Supervision of Health Care in accordance with applicable law.

Operational Documentation for Medical device
VITEK MS-CHCA – 411071 (01/2014)

Page 9 of 9

VITEK MS-CHCA matrix for use with VITEK[®] MS

IVD

The VITEK MS-CHCA matrix is used with the VITEK MS System to enable the identification of microorganisms. The matrix absorbs energy from the VITEK MS laser and transfers it to the microorganisms to enable ionization. The VITEK MS-CHCA matrix is used exclusively with the VITEK MS instrument as indicated in the VITEK MS Workflow User Manual ref. 4501-2233.

CONTENT OF THE KIT

Individual packaging:

REF 411071 VITEK MS-CHCA = 5 x 0.5 mL

COMPOSITION OF THE VITEK MS-CHCA MATRIX

Theoretical formula

Alpha-cyano-4-hydroxy-cinnamic acid	3.10 g
Ethanol (*)	25.57 g
Acetonitrile (*)	25.44 g
Solvent	to make 100 mL

(*) Medium F : HIGHLY FLAMMABLE

Medium Xn : HARMFUL

R11 : Highly flammable.

R20/21/22 : Harmful by inhalation, in contact with skin.

R36/37/38 : Irritating to eyes, respiratory system and skin.

S7 : Keep container tightly closed.

S16 : Keep away from sources of ignition - No smoking.

S26 : In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S36/37/39 : Wear suitable protective clothing, gloves and eye/face protection.

For more detailed information consult the material safety data sheet available on request.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use and microbiological control.
 - For professional use only.
 - For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
 - All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI[®] M29-A, Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
 - Do not use the VITEK MS-CHCA matrix after the expiration date.
 - Before use, check that the cardboard box is intact.
- Note: as this matrix is saturated, yellow crystals may be observed, but product performance is not affected.

STORAGE CONDITIONS

- Before use, store tubes at 2-8°C protected from light.
- The VITEK MS-CHCA matrix may be kept for up to 1 week after the tube has been opened.

QUALITY CONTROL

The VITEK MS-CHCA matrix is systematically quality controlled at various stages of its manufacture.

For users who wish to perform their own quality control tests, follow the test procedure described in the VITEK MS Workflow manual using the following strains:

- *Enterobacter aerogenes* ATCC[®] 13048[™]
- *Candida glabrata* ATCC[®] MYA-2950[™]
- *Escherichia coli* ATCC[®] 8739[™] (used as calibration strain)
- Negative test (matrix alone with no organism)

The criterion for quality control conformity is identification to species level with high level of confidence for the calibration and QC strains and no identification for the negative test.

WASTE DISPOSAL

The matrices must be discarded after use according to good laboratory practice.

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
LOT	Batch code
	Consult instructions for use
Rx only	For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

BIOMERIEUX, the blue logo and VITEK MS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

Any other name or trademark is the property of its respective owner.

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


Sample Preparation

4

IMPORTANT: Sample preparation is done in conjunction with the VITEK MS Prep Station. Review the VITEK MS Prep Station section of this manual (Section 2, Samples) before preparing the samples. The following section describes the process of placing samples on the VITEK MS-DS target slides in detail. Section 5, Preparation of VITEK MS-DS target slides and VITEK® 2 cards on the VITEK MS Prep Station software describes entry of sample information through the VITEK MS Prep Station.

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For measurement in the VITEK MS MALDI-TOF mass spectrometer, both the *Escherichia coli* ATCC® 8739 calibration strain and the sample must be freshly cultivated as described in the Sample section of this manual (See Section 2, Samples). The microorganisms to be identified must first be isolated on a suitable culture medium, according to established protocols. It is preferable to use colonies that are approximately three millimeters in diameter. However, it is also possible to use colonies that are considerably smaller or larger.

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Inactivation

CAUTION:
The inactivation step described below has to be performed using a microbiological safety cabinet.

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CAUTION:
When the target slide has been prepared it must be tested within 2 hours.

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Approved by:

Full name: Stéphanie HOMOLLE
Position: International Regulatory Affairs Specialist
Company Name: bioMérieux S.A.
Signature: *[Handwritten Signature]*
Date: *January, 31st, 2014*
Stamp:

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
Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Name of the device (IVD MD)	X	Yes	VITEK MS-FA Top of the first page of Package Insert (PI)
	Purpose of the IVD MD	X	Yes	Top of the first page of PI VITEK MS-FA is an auxiliary reagent for <i>in vitro</i> diagnostic and microbiological control. It is intended to be used for the preparation of samples for further analysis by means of MALDI-TOF (Matrix-Assisted Laser Desorption/Ionization – Time Of Flight), exclusively on VITEK MS mass spectrometer, manufactured by bioMérieux S.A., France (registration certificate No. FSZ 2012/12091 from 05.05.2012). A set of reagents (including VITEK MS-FA) is used to conduct a compulsory stage of preparation of microorganisms sample. The preparation of yeasts differs from the standard procedure in that VITEK MS-FA is applied to the sample prior to adding the matrix.
	Indication that the IVD MD is intended for use for clinical laboratory diagnostics	X	Yes	Section: Warnings and Precautions of PI
	Designated user	X	Yes	Section: Warnings and Precautions of PI

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Test principle		N/A	Not applicable in this case: VITEK MS-FA is an auxiliary reagent for <i>in vitro</i> diagnostic and microbiological control. It is intended to be used for the preparation of samples for further analysis by means of MALDI-TOF (Matrix-Assisted Laser Desorption/Ionization – Time Of Flight), exclusively on VITEK MS mass spectrometer, manufactured by bioMérieux S.A., France (registration certificate No. FSZ 2012/12091 from 05.05.2012). A set of reagents (including VITEK MS-FA) is used to conduct a compulsory stage of preparation of microorganisms sample. The preparation of yeasts differs from the standard procedure in that VITEK MS-FA is applied to the sample prior to adding the matrix. For more information, see: • Section 4 of User Manual - 4501-2233 - C - en - NA - VITEK MS Clinical Use – Workflow • Section 4 of User Manual - 4501-2233 - C - ru - NA - VITEK MS Clinical Use – Workflow
	Description of reagents, calibrators and control materials.		N/A	Section: Content of the kit of PI The kit does not contain any calibrators and control materials.

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Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Warnings and/or precautionary measures associated with potentially infectious materials contained in the IVD MD.	X	Yes	- Section: Warnings and Precautions of PI <i>Note: VITEK MS-FA does not contain any potentially infectious materials. However all specimens, microbial cultures and inoculated products to be analyzed on VITEK MS instrument should be considered infectious and handled appropriately.</i>
	As applicable: - special requirements for the premises - special training - special qualification of the expected user of the IVD MD.	X	Yes	Section: Warnings and Precautions of PI
	- Conditions for collection, processing and preparation of samples - Information about stability of the analyzed samples		N/A	Not applicable in this case: the VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method.
	Information required to verify whether the IVD MD is installed correctly and is ready to for safe operation	X	Yes	Section: Warnings and Precautions of PI
	Recommendations regarding the quality control procedures, where applicable.	X	Yes	Section: Quality Control of PI

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 1: Instructions for Use	Information about metrological traceability of the values assigned to the calibrators and the materials for control of correctness, including an indication to the applicable reference materials and/or reference measuring procedures of higher levels.		N/A	Not applicable in this case: The kit does not contain any calibrators and control materials.
	Analysis procedure, including calculations and interpretations of the results, and, where applicable, information about the feasibility of confirmation tests.		N/A	Not applicable in this case: - VITEK MS-FA does not have any measuring function. - VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method.
	Analytical parameters, such as sensitivity, specificity, correctness, repeatability, reproducibility, detection limits and measurement range, including information regarding control of known interfering agents, restrictions of the method and use of available reference materials and measurement methods.		N/A	Not applicable in this case: the VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method.
	Clinical efficacy parameters, such as diagnostic sensitivity and diagnostic specificity, as applicable.		N/A	Not applicable in this case: the VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method.
	Biological reference interval, where applicable.		N/A	Not applicable in this case: the VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method.

Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Information about interfering substances or restrictions associated with the sample, which can impact the study result.	N/A	Not applicable in this case: the VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method.
Warnings and/or precautionary measures associated with the safe disposal of the IVD MD and accessories, if any	X	Section: Waste disposal of PI
Name of the manufacturer and/or its authorized representative, as well as their addresses and contact details, from which technical advice and support can be obtained.	X	Bottom of last page of PI
Instruction of the user of the necessity to furnish a notice containing the respective information to the Federal Service on Surveillance in Healthcare of the Russian Federation according to the applicable laws in case of detection of any side effects which are not specified in the instruction for use or the operation manual of the medical device, the adverse reactions during its use, the peculiarities of interaction of medical devices one with another, the facts and circumstances posing threat to the life and health of citizens and healthcare professional in the process of application and operation of the medical devices.	X	See bottom of last page of Technical and Operational Documentation for VITEK MS-FA.

Content of Technical Documentation according to Appendix G of the "Methodical Recommendations Regarding the Procedure for an Expert Examination of Quality, Efficacy and Safety of Medical Devices"		Information required	Does it appear in the Package insert? (Yes/No)	Where to find the information?
Section 2: Additional instructions of use for the IVD MD intended for self-testing or testing near the sickbed	Detailed information about the analysis procedure, including preparation of reagents, taking/preparation of a sample, and information regarding the procedure for conducting the test and calculating the results.		N/A	Not applicable in this case: VITEK MS-FA matrix is not intended for self-testing nor for near patient testing.
	- Recommendations regarding the user's actions in case of a positive, negative or indefinite result; - Information about the test restrictions and the possibility of obtaining falsely positive or falsely negative results, as well as information regarding any factors which may influence the test result (e.g. age, sex, menstrual periods, infections, exercises, diets or medicine administration).		N/A	Not applicable in this case: VITEK MS-FA matrix is not intended for self-testing nor for near patient testing.
	For IVD MD intended for self-testing - an indication to the fact that the user should not make medical decisions without the prior consultation with the healthcare professional.		N/A	Not applicable in this case: VITEK MS-FA matrix is not intended for self-testing nor for near patient testing.

On issues related to the quality of medical products "VITEK MS-FA", and in case of detecting any side effects specified neither in the instruction for use nor the operating instructions on the medical devices, of adverse reactions when used, interaction of medical devices, as well as facts and circumstances endangering the life and health of citizens and medical professionals during the use and operation of medical devices please contact Authorized representative of the manufacturer:

LLC "bioMerieux RUSS"

REF 411072
VITEK MS-FA for
spectrometry is a meth
VITEK MS
spectroscopic mass. VITEK MS
Matrix-Assisted Laser
desorption/ionization (MALDI)
method (Matrix-Assisted Laser
desorption/ionization) is indi
cated in the VITEK MS
manual.
CONTENT OF THE KIT
Individual packaging: VITEK
REF 411072
COMPOSITION OF THE KIT
Theoretical formula
Formic acid (1)
Medium C: COF
R34 : Causes b
S26 : In case o
S27 : Take o
S28 : After r
and plenty
S36/37/39
eyeface
For more
sheet av
WARNING
• For i
• For
• For
de p

115114, Moscow
Derbenevskaya Street
Building 20
Structure 11

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Operational Documentation for Medical device
VITEK MS-FA – 411072 (01/2014)

Page 9 of 9

REF 411072

9300853D - en - 2013/07

EN

VITEK MS-FA for use with VITEK[®] MS

IVD

Mass spectrometry is a method for the detection and identification of molecules of interest by measuring their monoisotopic mass. VITEK MS-FA is used to pretreat yeasts to facilitate their identification according to the MALDI-TOF method (Matrix-Assisted Laser Desorption Ionisation-Time Of Flight). It is used exclusively with the VITEK MS instrument as indicated in the VITEK MS Workflow User Manual ref. 4501-2233.

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S27 : Take off immediately all contaminated clothing.

S28 : After contact with skin, wash immediately with soap and plenty of water.

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WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use and microbiological control.
- For professional use only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI[®] M29-A, Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use the VITEK MS-FA after the expiration date.
- Before use, check that the cardboard box is intact.

STORAGE CONDITIONS

- Before use, store the tubes at 2-8°C.
- VITEK MS-FA may be kept for up to 2 weeks after the tube has been opened.

QUALITY CONTROL

VITEK MS-FA is systematically quality controlled at various stages of its manufacture.

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- *Candida glabrata* ATCC[®] MYA-2950[™]
- *Escherichia coli* ATCC[®] 8739[™] (used as calibration strain)
- Negative test: VITEK MS-FA + VITEK MS-CHCA matrix (REF 411071) (matrix alone with no organism)

The criterion for quality control conformity is identification to species level with high level of confidence for the calibration strain and *Candida glabrata* strain and no identification for the negative test.

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INDEX OF SYMBOLS

Symbole	Signification
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
LOT	Batch code
	Consult instructions for use
Rx only	For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner

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BIOMERIEUX, the blue logo and VITEK MS are used, pending and/or registered trademarks belonging to bioMérieux or one of its subsidiaries, or one of its companies.

CLSI is a registered trademark belonging to Clinical Laboratory and Standards Institute, Inc.

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


Sample Preparation

4

IMPORTANT: Sample preparation is done in conjunction with the VITEK MS Prep Station. Review the VITEK MS Prep Station section of this manual (Section 2, Samples) before preparing the samples. The following section describes the process of placing samples on the VITEK MS-DS target slides in detail. Section 5, Preparation of VITEK MS-DS target slides and VITEK[®] 2 cards on the VITEK MS Prep Station software describes entry of sample information through the VITEK MS Prep Station.

This section explains how microbiological samples are prepared for MALDI-TOF mass spectrometry with the VITEK MS system.

For measurement in the VITEK MS MALDI-TOF mass spectrometer, both the *Escherichia coli* ATCC[®] 8739 calibration strain and the sample must be freshly cultivated as described in the Sample section of this manual (See Section 2, Samples). The microorganisms to be identified must first be isolated on a suitable culture medium, according to established protocols. It is preferable to use colonies that are approximately three millimeters in diameter. However, it is also possible to use colonies that are considerably smaller or larger.

Suitable colony size

18 to 24-hour culture

Note: The *E.coli* ATCC 8739 strain for the calibration should be incubated for 18 to 24 hours at 35°C +/- 2°C on blood agar.

Preparation of the VITEK MS-DS target slide

Note: Powder-free gloves should be used when handling VITEK MS-DS target slides.

Avoid touching the spots on the slide.


Handle the slide by the end of the slide.

Do not write or mark on the slide.

Before adding samples to the VITEK MS-DS target slide, place the VITEK MS-DS target slide on a flat surface, such as a silicone mat, that will prevent it from sliding.

The VITEK MS-DS target slide used for the VITEK MS is the disposable VITEK MS-DS target slide with 48 sample positions. This VITEK MS-DS target slide is for single-use only. The sample positions are divided into three sections. Each section is called an acquisition group. The three small spots in the center of each acquisition group must not be used for identification of microorganisms. These three small spots (xA1, xB1, xC1) are used for calibration.

Acquisition Groups

Position of *E. coli* ATCC® 8739 calibration strain on the VITEK MS-DS target slide

The *E. coli* ATCC 8739 strain has to be deposited on positions: xA1, xB1, xC1, depending on the number of samples tested (1 calibrator per acquisition group of 16 spots).

Follow the directions for Preparation of Bacteria to prepare the *E. coli* ATCC 8739 calibration strain.

Preparation of bacteria

Note: Powder-free gloves should be used when handling VITEK MS-DS target slides.

Using a 1 µl loop, pick up part of a suitable colony.

Make sure not to pick up any agar with the colony as this would have a negative impact on the spectrum.

Some microorganisms, such as Streptococci, grow in very small colonies. In such cases, pick up several similar colonies and deposit them on one spot.

Note: Do not deposit more than 1 (Single Deposit Mode) or 2 (Duplicate Deposit Mode) samples at a time before adding matrix to the target spots.

Apply the sample to the center of the spot. Take special care not to apply too much. The appropriate amount on a spot is shown in the picture below. If not enough colony was deposited, take a new loop and add another part of the colony.

Too much Too little OK OK

Smear a thin layer of the sample on the spot, then discard the used loop.

Open the microtube containing the VITEK MS-CHCA matrix. Add precisely 1.0 µl of matrix to the center of the spot. Discard the pipette tip. The matrix/microorganism suspension is then left to dry completely.

It is strongly recommended to deposit the matrix as indicated in the VITEK MS Prep Station workflow for bacteria:

- 2 by 2 in duplicate mode,
- 1 by 1 in single mode.

Note: The VITEK MS-CHCA matrix contains a high concentration of organic solvents; it is recommended to close the tube after dispensing to avoid evaporation.

Proceed likewise for a second sample...etc...

After about five minutes, check for crystal formation on the spots. Inexperienced users may find using a magnifying glass helpful. Measurement will succeed only if matrix crystals have become visible as a yellowish film. Ideally, most of the spot's surface will be coated with crystals.


When the VITEK MS-DS target slides have been prepared, they must be tested within 48 hours.

Prepared slides should be stored at room temperature until they are tested.

Preparation of yeasts

Note: Powder-free gloves should be used when handling VITEK MS-DS target slides.

The preparation of yeasts differs from the standard procedure in that VITEK MS-FA is applied to the sample prior to adding the VITEK MS-CHCA matrix.

Yeast immediately after VITEK MS-FA treatment

Begin with applying the sample to the appropriate VITEK MS-DS target slide position as in the standard procedure. Smear a thin layer of colony. Add 0.5 μ l of VITEK MS-FA. Deposit a second sample and using a fresh pipette tip, add 0.5 μ l VITEK MS-FA ...etc...

Note: The VITEK MS-CHCA matrix contains a high concentration of organic solvents; it is recommended to close the tube after dispensing to avoid evaporation.

For efficient extraction, it is important to allow for evaporation of the VITEK MS-FA before adding the matrix. Depending on the air humidity and other factors, this usually takes between one and three minutes.

Add 1.0 μ l of VITEK MS-CHCA matrix following the same procedure as for bacteria.

On a single VITEK MS-DS target slide, it is possible to prepare samples using both procedures (bacteria and yeasts).

CAUTION:

Do not treat the *E.coli* calibration strains on positions xA1, xB1, xC1 with VITEK MS-FA.

When the VITEK MS-DS target slides have been prepared, they must be tested within 48 hours.

Prepared slides should be stored at room temperature until they are tested.

Preparation of Mycobacteria

Note: Powder-free gloves should be used when handling VITEK MS-DS target slides.

The preparation of Mycobacteria differs from the standard procedure in that an inactivation step is necessary prior to adding the VITEK MS-CHCA matrix.

Inactivation

CAUTION:

The inactivation step described below has to be performed using a microbiological safety cabinet.

The microorganisms to be identified must first be isolated on a suitable culture medium according to established protocols.

For each sample to be identified, prepare a microtube (2 ml) containing 50 µl of Trifluoroacetic acid (TFA).

From the solid medium, take a minimum quantity of bacterial cells equivalent to 2/3 of a 10 µl calibrated loop and add directly into the TFA. In order to homogenize, the bacterial suspension has to be vortexed a few seconds. Then the inactivation step is started, wait 30 minutes until completion.

Suspension preparation

This inactivated suspension has to be diluted 1/10 by adding 450 µl demineralized and deionized water to the microtube. Vortex a few seconds. This suspension is now ready to be deposited on the VITEK MS-DS target slide.

Deposit

Deposit 1 µl of the prepared suspension on the appropriate spot using a micropipette. Using a fresh pipette tip, deposit another sample on the next spot ... etc. Allow to dry (5 to 10 minutes) and then add 1 µl of VITEK MS-CHCA matrix. Depending on the air humidity and other factors, this usually takes several minutes. It is important to add the matrix when the inactivated bacterial suspension has become completely dry. This is usually after one to three minutes.

CAUTION:

When the target slide has been prepared it must be tested within 2 hours.

Preparation of Moulds

CAUTION:

Use a safety cabinet, wear a laboratory coat, gloves and oversleeves.
Use the appropriate procedure / product to clean the material.
At the current time, only the Sabouraud Dextrose Agar medium has been validated.

- 1) Wet a swab with suspension medium and take a piece of moulds (spores or hyphes can be picked up).

Note: Do not remove agar.

- 2) Resuspend moulds in a 2 ml Eppendorf tube containing 300 µl API suspension medium.

- 3) Add 0.9 ml 99.8%, 96% or Absolute Ethanol.

- 4) Centrifuge for 2 minutes at 10,000 rpm.

Note: Please note that depending on the instrument, the centrifugation time has to be adapted (for example 2:30 min in order to have 2 min at 10.000 rpm). Moreover, if you need to take the tubes out of the safety cabinet for the centrifugation steps, carefully clean the tubes with sporicid product in order to avoid contamination.

- 5) Discard the supernatant using a pipette.

- 6) Add 40 μ l Formic acid 70%* and mix (vortex).

Note: *Formic acid 70% must be prepared each time an identification of moulds is required. The preparation must be done under a chemical safety cabinet as follows: 3 volumes of Suspension Medium + 7 volumes of pure formic acid.

- 7) Add 40 μ l Acetonitril and mix (vortex).

- 8) Centrifuge for 2 minutes at 10.000 rpm.

- 9) Deposit 1 μ l supernatant.

- 10) Let dry.

- 11) Add 1 μ l CHCA matrix.

- 12) Let dry.

CAUTION:

When the target slide has been prepared it must be tested within 2 hours.

Sample preparation failures

If none of the VITEK MS-DS target slide position(s) for a sample yields satisfactory identification results, check the appearance of the deposits (the deposit has been spread outside the well limits, the deposit is too dry and peels off the VITEK MS-DS target slide etc.).

If contamination or poor crystal formation could have been the cause of the unsatisfactory results or if no matrix was added, take a new VITEK MS-DS target slide and repeat the procedure.

Contamination

The VITEK MS-DS target slides must be stored in a dust-free environment and are designed for single-use only.

Powder-free gloves should be used when handling VITEK MS-DS target slides.

Do not touch the upper face of the VITEK MS-DS target slide before use.

Take care not to contaminate the VITEK MS-DS target slide with bacteria and avoid liquid spills.

Do not use loops that may have come into contact with anything other than the sample strain to be tested.

Picking up agar along with the colony may lead to poor identification results.

Be careful that only fresh, sterile pipette tips come into contact with the matrix or the VITEK MS-FA.

Inappropriate quantity of colony

Automated measurement generally yields high quality mass spectra for a wide range of cell material amounts applied to the VITEK MS-DS target slide positions. Nevertheless, an insufficient or excessive quantity of colony may result in failed spectra or in spectra of low quality. Applying an insufficient quantity of colony usually results in no spectra being acquired. Applying too much colony is the common cause of suboptimal performance of the system. If an excessive quantity of colony is applied, the cells may not suspend well in the matrix suspension. This hampers the extraction process and the subsequent crystal formation by the matrix.

Note that colonies on agar plates may vary considerably with respect to size, texture or viscosity, requiring adaptation of the technique for application.

Inappropriate application of the matrix

Make sure to always add the appropriate volume of matrix. When using VITEK MS-CHCA matrix, 1 μ l is added to each spot, whereas only 0.5 μ l is added when using VITEK MS-FA for yeast extraction. If by mistake no matrix is added to a spot, no identification result will be obtained.

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