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Dear Customer,

Z4L medolight is manufactured with the latest production methods and quality components to ensure trouble-free performance. If your appliance still needs warranty service, please send it together with this guarantee certificate to your Z4L medolight dealer from whom you bought it or to authorized Zepter service.

The Guarantee

With this Consumer Guarantee, ZEPTER International Poland Sp. z o.o. guarantees this product to be free of defects in materials and workmanship at the time of its original purchase or at the time it was taken on hire purchase terms by the consumer from the retailer for the period of TWO (2) years. If during this period of guarantee the Z4L medolight Therapy Device proves defective due to improper materials or workmanship, the Retailer, Authorized Service Agents or Authorized Servicing Dealers will, without charge for labour or parts, repair or replace this product or its defective parts on the conditions explained hereafter.

Conditions

- I. The guarantee is valid when presented together with the original invoice/receipt issued to the consumer by the retailer/sales rep, and if the guarantee card states a) the purchaser's name, b) the retailer's name and address, c) the model name and serial number of the purchased product and d) the date of purchase of the product. ZEPTER International Poland Sp. z o.o. reserves the right to refuse guarantee service if this information has been removed or changed after the original purchase of the product by the consumer from the retailer.
- II. This product will not be considered defective in materials or workmanship if it has been adapted, changed or adjusted in order to conform to the national or local technical or safety standards in any country other than the one for which the product was originally designed and manufactured for. This guarantee does not cover reimbursement for a) such adaption, or adjustments or attempts to do so, whether properly performed or not, or b) any damage resulting thereof.
- III. This guarantee covers none of the following:
 - A. Periodic check-up, maintenance, and repair or replacements of parts due to normal wear and tear;
 - B. Home service transport costs and other costs and risks of transport relating directly or indirectly to the guarantee of this product.
 - C. Damage to this product resulting from:
 1. Abuse and misuse, including but not limited to:
 - a) the use of this product for other than its intended purposes or not in accordance with manufacturer's instructions for the proper use and maintenance of this product, and OR
 - b) the installation of this product in a manner inconsistent with the technical or safety standards in the country where the product is used;
 2. Any repairs carried out by unauthorized Persons or Service Stations
 3. Accidents, acts of God or any cause beyond the control of ZEPTER International Poland Sp. z o.o., including but not limited to lightning, water, fire, public disturbances and improper ventilation.
- IV. This guarantee only covers Z4L medolight units which were imported and sold by an authorized Z4L medolight retailer/sales representative.
- V. The guarantee claim must be made via your retailer/ send directly to authorized service centre.
- VI. This guarantee does not affect the consumer's statutory rights under applicable national legislation in force, nor the consumer's rights against the retailer arising from the sales/purchase contract.

In the absence of applicable national legislation, this guarantee will be the consumer's sole and exclusive legal remedy, and neither ZEPTER International Poland Sp. z o.o. nor its Subsidiaries or Distributors shall be liable for any incidental or consequential damages for breach of any expressed or implied guarantee of this product.

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