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SINGAPORE ACADEMY OF LAW

I, Low Hui Min, Chief Financial Officer, Singapore Academy
of Law, Republic of Singapore, hereby certify that Victor Koh Koon
Hor is a duly appointed Notary Public Practising in Singapore, and
that the signature appearing at the foot of the annexed Notarial
Certificate dated 17th February 2016 is the signature of the said
Victor Koh Koon Hor.

Dated at Singapore this 17th day of February 2016.

LOW HUI MIN
CHIEF FINANCIAL OFFICER
SINGAPORE ACADEMY OF LAW

16021261

Certified true signature

EJILANE

18 FEB 2016

TO ALL TO WHOM THESE PRESENTS SHALL COME

I, **VICTOR KOH KOON HOR**, Notary Public, duly authorised admitted appointed and practising in the Republic of Singapore, DO HEREBY CERTIFY that the document hereunto annexed is a true copy of the original document of which it purports to be a copy, I having carefully collated and compared the true copy with the original and found the true copy to agree with the original.

IN TESTIMONY WHEREOF I, the said Notary have hereunto subscribed my name and affixed my seal of office this 17th day of February, 2016.

WHICH I ATTEST

Victor Koh Koon Hor


NOTARY PUBLIC,
SINGAPORE

Infrared Thermometer MediCare MC 302

Instruction Manual


MED**CARE**

Approved Signature: _____

Name: Amos Sim

Position: Business director,

Management Representative

CERTIFIED TRUE COPY

Victor Koh Koon Hor

17 FEB 2016

Contents

1. Product's presentation and its classification	
2. Recommendations for use, indications, contraindications	
3. Principle of operation	
4. Characteristics	
5. Technical characteristics	
6. Device description	
7. Component parts	
8. Buttons functional description	
9. Display message description	
10. Explanation of symbols used on the consumer package and label	
11.1 Battery installation/replacement	
11.2 Setting of units of temperature measurement - F-1	
11.3 Setting of a testing mode - F-2	
11.4 Setting of audio parameters- F-3	
11.5 Setting of temperature threshold values - F-4	
11.6 Setting of a general deviation value - F-5	
11.7 Laser indicator switch on - F-6	
11.8 Exit	
12. Recommendations for use	
12.1 Basic instructions	
12.1.1 Body temperature measurement	
12.1.2 Object temperature measurement	
12.2 View results	
12.3 Environment temperature measurement	
12.4 Measurement in the "scan function"	
12.5 Parameters beyond the measurement range	
13. Safety information	
14. Care and cleaning	
15. Troubleshooting	
16. Conditions of storage, transportation and disposal regulations	
17. Environment protection requirements	
18. Quality commitment and post-marketing maintenance	
19. Compliance with standards	
20. Accuracy check	
21. Manufacturer	
22 Warranty card	

The instruction manual contains information on the operation of the infrared thermometer MediCare MC 302. Read this instruction before using the thermometer.

1. Product's presentation and its classification

Infrared Thermometer Medicare MC 302 is designed to measure the surface temperature of the body or an object by a non-contact manner. Thermometer Medicare MC 302 is an infrared thermometer, the principle of operation is based on the measurement of infrared radiation. The thermometer is characterized by the convenient use, robust design, high-precision measurement and interference protection. For fast and accurate measurement of the body temperature in one moment, you just need to direct the thermometer on the forehead area or an object, the temperature of which is planned to be measured and press the button. It is ideal for measuring body temperature in children, even during their sleep. You can also measure the room temperature, the temperature of various surfaces, bathing water, baby food, etc. Thermometer is widely used domiciliary, in schools, hospitals, when testing food, in agriculture and industry.

2. Recommendations for use, indications, contraindications

Infrared Thermometer Medicare MC 302 is designed for non-contact body temperature measurement of adults and children («Body temp» function), it can be used for an object temperature measurement or room temperature measurement («Object temp» function).

The infrared thermometer is not designed for temperature measurement by rectal, oral or axillary (underarm) methods.

The following factors influence the temperature of a healthy person:

- Specific features of metabolism.
- Age (body temperature of infants and young children is higher and it decreases with age. The temperature fluctuations in children occur faster and more often).
- Clothes worn by a person.
- Ambient temperature.
- The time of a day (body temperature is lower in the morning and it decreases by the end of a day).
- Prior physical exercise.
- Method of measurement.
- The phase of the menstrual cycle.

Indications: for fast and accurate body temperature measurement.

Contraindications: there are no known contraindications.

Possible side effects are not identified.

Note: we recommend you to know your individual normal body temperature when you feel good, in order to interpret the data properly which differ from the average data (e.g. fever). Depending on the type and thickness of the skin there may be a difference in the body temperature readings.

Warning:

in case your body temperature is above 38 degrees, we recommend you immediately contact your physician.

The temperature measurement depends on the body part, the temperature of which is measured. In healthy people, such fluctuations for different body parts can range from 0.2 °C to 1 °C. You can't compare the temperature measured by various thermometers. When measuring the temperature, take into account a thermometer you use and in what body part you measure the temperature, and always inform your doctor about it and take it into account when self-testing. The ambient temperature also affects the skin temperature.

Keep in mind that the infrared thermometer Medicare MC 302 and a person or an object the temperature of which will be measured, must be present in a room, where the temperature will be measured, for at least 10 minutes before the measurement.

For achieving correct and constant results it is recommended to check and, if necessary, set the device when using the thermometer for the first time (clause 9.5 of this instruction).

3. Principle of operation

Any object at above zero centigrade generates a certain proportion of infrared radiant energy depending on its own temperature. The emission energy distribution wavelength depends on the surface temperature. The body temperature (35~43 °C) can release an infrared ray at the range of 9~13 μm. The principle of operation of the infrared thermometer Medicare MC 302 is based on the measurement of the absolute value of the electromagnetic emission amplitude from the object in the infrared spectrum and the subsequent conversion of the measured values in the temperature. The infrared ray focused by an optical system is incident on the sensor (primary pyrometrical transducer) forming an electrical output signal proportional to the temperature value of the measurement object. This signal passes through an electronic transducer (secondary pyrometrical transducer), enters the measuring counting unit and is processed in it. The result is displayed.

4. Characteristics

- Easy to use, allows you to get an accurate and fast result.
- It is designed to test forehead temperature of a person but it can be also adjusted to test environmental temperature and object surface temperature.
- The use of high-precision infrared sensor provides reliable performance.
- Fast voice message concerning a very high temperature (a value can be set in advance).
- It saves automatically up to 64 measured values.
- A large backlit LCD screen.
- Availability of Fahrenheit (°F) and Celsius (°C) temperature scales.
- The object temperature mode can measure the temperature of an object within the range of 0-118 °C when the emissivity is 0.95.
- Lock automatic scan function.
- Automatic switch off function to extend the battery service life.

5. Technical characteristics


Measurement method	Non-contact infrared measurement
Measurement units	° Celsius (°C) or ° Fahrenheit (°F)
Operation environment	Environment temperature: 10°C~40°C
	Relative humidity: от 10% до 80%
	Atmospheric pressure: от 900 до 1100 gPa
Size	225 mm x 54 mm x 78 mm (Length x Width x Height)
Weight	175 g (without batteries)
Display resolution	0.1 °C/ 0.1 °F
Temperature measurement range	Body temperature: 32.0°C~42.9°C (89.6°F~109.2°F)
	Object surface temperature: 0°C~118°C (32°F~244.4°F)
Measurement error	Under normal environment: ±0.2°C
Power consumption	Voice message: ≤50 mA
	No voice message: ≤15 mA
	Sleep mode: ≤10 mA
Measurement time	≤0.8 sec
Testing distance	5 cm~8 cm from the point of measurement
Emissivity	0.95
Automatic switch-off time	10 sec
Memory function of body temperature measurements	Automatic save of the last 64 measured values
Power supply	DC3V(two batteries 1.5 B type AAA)
Laser safety class	2
Contact type	B
Software	XZ-THIR-002 class A.dated 2011-06-18
Service life	3 years

Date of manufacture of the device is indicated in the serial number.

The manufacture year is coded in a letter code A = 1, B = 2, C = 3, D = 4, E = 5, the first and second digits indicate the month of the manufacture, the third and fourth digits indicate the manufacture date. The last four digits indicate the serial number of the device.

Example: if serial number is B0AD04060001, it means that the manufacture date: year 2014, month 04, date 06.
Serial number of a device is 0001.

6. Device description


The materials of which the infrared thermometer is made

Part	Material
Battery lid, lower casing, switch, buttons, front cover	ABS plastic
Upper casing, base	ABS plastic
Lens, display	Polycarbonate

7. Component parts


- Infrared Thermometer Medicare MC 302
- Base
- Instruction manual
- Battery AAA (1,5 V) – 2 pcs.

8. Buttons functional description

Button	Functional description
+	Used when setting the parameters, increasing the parameter or selecting the relevant parameter.
	Allows performing an upward page turning to check the saved values.
—	Used when setting the parameters, decreasing the parameter or selecting the relevant parameter.
	Allows performing a downward page turning to check the saved values.
SET	Used when setting the parameters. To enter the set up mode, press and hold the SET button for three seconds.
	Allows checking of the previous measurements. To check the previous records press the SET button.
SCAN	Used to turn on "Scan functions" when making measurements. To exit the "Scan functions" press the switch for one second.


9. Display message description


Sign	Function description
	Battery power level. When it appears or flashing, the battery is discharged.
	Laser indicator status
Body temp	Body temperature measurement mode. The displayed readings are body temperature values.
Object temp	Object surface temperature measurement mode. The displayed readings are object surface temperature values.

°C	Temperature measurement units in Celsius (°C).
°F	Temperature measurement units in Fahrenheit (°F).
	Voice message
M	Entrance of the required status to check past records.
	Measurement performance
Scan	«Scan function» is switched on during the measurement

10. Explanation of symbols used on the consumer package and label

Symbol	Explanation
CE 0499	Compliance with the directive CE 0499
	Contact detail of type B
	Caution! Read the instruction carefully before use.

	Special disposal. Do not dispose of with household waste. To avoid doing harm to the environment, it is necessary to separate the object from the normal waste and dispose of it in the safest way, for example, put in a special place for disposal.
	Laser emission
	Manufacturer
	Authorized Representative in the European Community

11. Preparation for the operation

After unpacking the device, insert the batteries included in the kit before the first use. The battery compartment is located on the bottom of the device. To set the thermometer's parameters, use the «SET» button.

11.1 Battery installation/replacement

1. Remove the battery cover.

2. Insert four batteries "AAA" according to the diagram on the compartment and place the cover back.

When the battery charge is low, the sign  appears on the display. In this case, it is necessary to replace the batteries as soon as possible. When the flashing symbol appears on the display  even after the battery replacement, it is necessary to make sure of the accuracy of measurement.

To replace batteries, open the cover of the battery compartment. Remove the used batteries. Insert two new alkaline batteries AAA, observing the correct polarity, and then close the cover. If the thermometer is not used for a long period of time - store it without batteries.

Batteries must be disposed of in accordance with state regulations for the disposal of batteries and supply elements.

11.2 Setting of units of temperature measurement - F-1

- Press and hold the «SET» button for 2 seconds, F-1 appears on the display.
- Press "+" button to choose Fahrenheit °F as the measurement unit scale, 'F' sign appears on the display.
- Press "-" button to choose Celsius °C as the measurement unit scale, 'C' sign appears on the display.
- When there is a necessity to enter setup F-2 mode, press the «SET» button.

11.3 Setting of a testing mode - F-2

- Press and hold the «SET» button for 2 seconds, F-1 appears on the display.
- When pressing the «SET» button for the second time, F-2 appears on the display.
- To select the object surface temperature measurement mode, press the "+" button, the message 'object temp' appears on the display.
- To enter the body temperature measurement mode, press the "-" button, the message 'body temp' appears on the display.
- When there is a necessity to enter setup F-3 mode, press the «SET» button.

Caution:

To obtain accurate results in the 'body temp' mode, it is necessary to direct the device precisely on the forehead area, previously having removed all possible obstacles such as hair, sweat and others.

11.4 Setting of audio parameters- F-3

- Press and hold the «SET» button for 2 seconds, F-1 appears on the display.
- Press the «SET» button until F-3 sign appears on the display.
- To switch off the audio signals, press the " - " button until 'CLOS' sign appears on the display.
- To switch on the audio signals, press the "+" button until 'BEEP' sign appears on the display.

11.5 Setting of temperature threshold values - F-4

- Press and hold the «SET» button for 2 seconds, F-1 appears on the display.
- Press the «SET» button until F-4 sign appears on the display;
- To increase the threshold value by 0.1 °C, press the "+" button once (the maximum value is 42 °C);

- To decrease the threshold value by 0.1°C , press the “-” button once (the minimum value is 37°C);
- When there is a necessity to enter setup F-5 mode, press the «SET» button.

If the temperature exceeds the threshold value, the audio signal appears which indicates that the temperature has exceeded the threshold limits. By default, the threshold temperature is 38°C .

Press and hold the button ‘+’ or ‘-’ for quick adjustment of values.

11.6 Setting of a general deviation value - F-5

- Press and hold the «SET» button for 2 seconds, F-1 appears on the display.
- Press the «SET» button till F-5 sign appears on the display.
- Press the button “+” to increase the deviation value by 0.1°C (maximum 3.0°C);
- Press the button “-” to decrease the deviation value by 0.1°C (minimum 0.0°C);
- When there is a necessity to enter setup F-6 mode, press the «SET» button.

Before use it is necessary to adjust the acceptable values of the actual deviation of the measured temperature in accordance with an interval of study and the environmental conditions from 0.0°C to $+3.0^{\circ}\text{C}$. The factory setting is the value of 0.8°C .

Example: measure the body temperature using a mercury thermometer. Let's suppose that you have got the result of 36.5°C . Measure the temperature with the infrared thermometer, keeping it at a distance of 5-8 cm from the forehead of the patient. Remove any possible obstacles from the patient's forehead (hair, sweat, etc.). If when measuring you've got a temperature of 36.50°C , the thermometer operates properly. If the temperature is different, e.g. 35.40°C (the difference is 1.10°C), it is necessary to set the thermometer. Using the setting of a general deviation value «F-5», adjust the deviation up to 1.10°C . After the setting, measure the temperature with a non-contact thermometer again to verify the correctness of the measurement.

11.7 Laser indicator switch on - F-6

- After pressing the «SET» button for 2 seconds, F-1 appears on the display;
- Press the «SET» button till F-6 sign appears on the display;
- Press the button “-” to turn off the laser beam indicator;
- Press the button “+” to turn on the laser beam indicator.

11.8 Exit

- To exit according to the setup parameters, press the button SET continuously until "SAVE" sign appears.
- Exit is carried out automatically if you do not press any button within twelve seconds.

Caution:

Only when pressing the SET button until the appearance of the sign «SAVE», all the set parameters can be saved, otherwise the operations described in clauses 9.1-9.6 will be invalid.

12. Recommendations for use

12.1 Basic instructions

12.1.1 Body temperature measurement

- Direct the laser indicator of the device on the central part of the forehead area at a distance of 5-8 cm;
- Press the switch to start the measurements;
- When measuring, the sign  appears on the right bottom corner of the screen;
- Release the switch and the measurement result will appear on the display.

The measured value will be saved automatically. See the picture below:

Body temp

36.8°C


- If within ten seconds there is no operation performed, the device will turn off automatically.

Caution:

- To ensure accurate test results, do not use a thermometer for ten minutes after replacing the battery (batteries) or when moving the device to a new environment or another room.
- Make sure before using that the measurement in the forehead area is not interfered with hair, sweat, hats and others.
- Deviations occur when measurements are not performed at a proper distance, or when there is a deviation from the central area of the forehead area of the patient.

12.1.2 Object temperature measurement

- Direct the laser indicator of the device on the tested surface at a distance of 5-8 cm to the measured object
- Press the switch to start the measurements;
- When measuring, the sign  appears on the right bottom corner of the screen;
- Release the switch and the measurement result will appear on the display. The measured value will be saved automatically. See the picture below:


- If within ten seconds there are no operation performed, the device will turn off automatically.

12.2 View results

- To view the stored measurements, press the SET button one time, the sign 'M' and the corresponding number of measurement will appear at the top right part of the display.
- To view the saved measurements, press '+' or '-' button. Every time you press '-' button, the next, earlier measured value will be displayed. When pressing "+" button, the following measured value will be displayed.

Display of
previous records
of body temperature

Request of previous records
measurement results

Serial number of a record.

Caution:

All test results are stored automatically after a measurement. The maximum number of stored records is 64. In the case when there are more than 64 records, a new value will automatically replace the previous one.

All previous records are marked with a serial number, for example, M 1 is the latest record, while M 64 is the first record.

12.3 Environment temperature measurement


To measure the ambient temperature, it is necessary to take thermometer off into an «Object temp» mode and not directing at any object, press the switch. The ambient temperature will be displayed.

12.4 Measurement in the "scan function"

See the procedure below:

a) Press the switch to start the thermometer, then press the «Scan» button to run the "Scan function" mode, the

sign «Scan» and  will be displayed, see the picture below:

b) To make a measurement direct the thermometer on the test surface or the forehead area of a person. To make a remeasurement it is not necessary to press the switch;

c) To stop «Scan function» press the switch and the «Scan» sign will disappear.

The Scan function is particularly useful in measuring the temperature of several objects or surfaces consistently.

Caution:

Please, stop the Scan function after the measurement to avoid a rapid battery discharge.

12.5 Parameters beyond the measurement range

- Body temperature measurement mode:

The sign «Lo» will appear on the display, if measurement results are below 32°C.

The sign «Hi» will appear on the display, if measurement results are above 42.9°C.

- Object temperature measurement mode:

The sign «Lo» will appear on the display «Lo», if measurement results are below 0°C.

The sign «Hi» will appear on the display, if measurement results are above 118°C.

Attention!

When the ambient temperature is less than 10.0°C or higher than 40.0°C, the wrong values may be displayed. At these temperatures, the measurement accuracy is not guaranteed.

13. Safety information

- The device must be used according to its intended purpose in accordance with this instruction.
- The device must be placed for at least 10 minutes in the room where the measurement will be performed before measuring.
- The device should be used according to the purposes described in this instruction manual.
- Using the device by children is prohibited.
- Keep out of children's reach.
- Check before every use if the lens of the thermal sensor is not damaged.
- The thermometer contains a laser indicator, never point a laser directly at eyes.
- Do not use optical instruments like telescope, etc. for investigating the laser beam.
- The device has been designed for a practical use, but it can't replace a visit to a doctor.
- Do not expose to the mechanical impact and do not use the device in case of its damage.
- Do not expose the thermometer to electricity.
- Do not expose to high and low temperatures ($> 50^{\circ}\text{C}$ or $< 0^{\circ}\text{C}$)
- Do not immerse the device in water or other liquids
- Do not use the device near sources of high electromagnetic emission, such as sources of wireless communications or cell phones.
- Keep away from water and thermal sources, including direct sunlight.
- Protect the device from extreme temperatures, humidity, dust and direct sunlight.
- Thermometer's probe is the weakest part of the thermometer. it must be handled carefully.
- Do not try to charge non-rechargeable batteries.
- In case of problems in the thermometer operation, contact the service center. Do not try to fix this problem by yourself.

14. Care and cleaning

- In order to ensure the accuracy of measurement and to avoid microbial contamination, every time before use and/or after use, wipe the thermometer lens with a soft cotton cloth or a cotton bud moistened with water or alcohol. Leave the device for 30 minutes after wiping until it gets dry.
- The casing of the thermometer should be wiped with a soft, dry cloth.

- Keep away from water and other liquids in a dry location away from dust and direct sunlight.
- Do not use abrasive cleansers and do not immerse the thermometer in water or other liquids.

15. Troubleshooting

Trouble description	Troubleshooting method
The measurement results are too low or too high.	1. Check if the deviation value is changed, according to the factory setting it is 0.8°C. 2. The thermometer can have errors due to environmental changes. If the environmental temperature changes considerably or if the thermometer is used firstly for high-temperature object and then for a very low-temperature one, the measurement difference will happen. The thermometer should be kept in a relatively stable environment during 10 to 30 minutes to reach a balance before starting a new measurement. Measurement distance is from 5 to 8 cm.
No responses after pressing a key	Do not dismantle and do not repair the thermometer by yourself. Contact the organization authorized to accept claims.
No display or values are distorted	Do not dismantle and do not repair the thermometer by yourself.
No voice signal	Enter F3 to check switch on/switch off parameters.
No response after pressing a control key	1. Batteries don't have enough power, and battery replacement is necessary. 2. Check if the battery polarity is placed in correct position.

If, in spite of the above mentioned recommendations, you can't achieve the correct results, stop the device operation and contact the company authorized to accept claims for a technical support.

Never try to repair the device by yourself. Any attempt to open the casing of the device and fix it, removes any warranty liabilities under the warranty repair of the device.

16. Conditions of storage, transportation and disposal regulations

The packed products should be stored in a room that has good ventilation system. Drops, strong pressings, precipitation and step movement should be excluded during the storage and transportation. It is necessary to indicate the direction of the location on the package.

Storage and transportation conditions

Environmental parameter	Value
Environmental temperature	From +10°C to +55°C
Relative humidity	from 10 to 85 % (without gas corrosion)
Atmosphere pressure	from 900 to 1100 GPa

The disposal of this device and used batteries should be carried out in accordance with national regulations for disposal of electronic products.

17. Environment protection requirements

The products of L-Tac Medicare Pte Ltd. company are designed and manufactured in accordance with applicable guidelines for the protection of the environment. Provided proper disposal, you can avoid the risk of environmental pollution.

18. Quality commitment and post-marketing maintenance

The manufacturer guarantees that the product complies with the technical specifications and performs its functions providing a correct use and observance of storage conditions and transportation conditions specified in the instruction manual.

The manufacturer provides warranty liabilities for 12 months from the date of purchase.

The warranty is valid only if there is a sales receipt or the warranty card with all information filled in with a seal of a commercial enterprise or a technical service center.

The guarantee is not performed in the presence of traces of mechanical impact, dents, cracks, chips, etc on the device (display included); in the presence of traces of opening the case, traces of attempts of repair in an authorized service center; in the presence of traces of ingress of moisture inside the casing, or exposure of the device to aggressive environment; in other cases of violation by a consumer of storage regulations, transportation and technical operation of the device, provided by the regulations, indicated in the manual.

The guarantee does not cover batteries and packaging.

If, during the warranty period, products are used not in accordance with the intended use or with violations of the instruction manual, repair is made by the customer.

The company authorized to accept claims:

Formula-FR LLC,
105043, Russia, Moscow, Pervomayskaya street, 35/18.
Tel. 8-926-208-71-61

Please, keep the warranty card and the purchase receipt for a future maintenance.

19. Compliance with standards

The high quality of the device is confirmed by documents.

In Russia

Registration certificate _____

Declaration of conformity.

Certificate of type of measuring devices of the Federal Agency on Technical Regulating and Metrology

In the European Union

The manufacture of devices is certified according to the international standard ISO 13485:2003.

The device conforms to the European Directive MDD 93/42/EEC, international standards EN 12470-5:2003, ISO 14971:2009, ISO 14151-1:2009, EN ISO 14151-2:2009, EN 980:2008, EN 1041-2008, IEC/EN 601-1-2:2007, IEC/EN60601-1:2005(3rd Ed) EN ISO 10993-1:2009, EN ISO 10993-5:2009, EN ISO 10993-10:2009.

20. Accuracy check

We recommend checking the accuracy every 2 years or after mechanical impact (e.g. drops).

The verification interval is 2 years.

More information about verification can be obtained in the local branch office of Rostest.

21. Manufacturer

«L-Tac MediCare Pte Ltd», Singapore

Manufacture site address

25 Mandai Estate, #07-11

Innovation Place Tower 1

Singapore 729930

Tel : +65 6863 8191

Fax : +65 6863 8190

www.medicare-l-tac.com.sg

Customer Service: +65 6334 2273

22 Warranty card

Serial number _____

Date of purchase _____

Name and address of the trade organization _____

Seller's signature _____ (stamp)

The product is checked and has no damages

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