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EC Declaration of Conformity

Asclepion Laser Technologies GmbH

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07747 Jena

Germany

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We declare the compliance of the medical device

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Essential requirements of Directive 93/42/EEC, Annex I, with all relevant harmonized standards

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with the requirements of the Medical Devices 93/42/EEC Annex II (registration no. 1901839-006-000 of the notified body "LGA – InterCert GmbH", Tillystr. 2, 90431 Nürnberg, no. 1275). ¹⁾ Any modification to the product, not authorized by Asclepion Laser Technologies GmbH, will invalidate this declaration.

Jena, June 28th 2013

Asclepion Laser Technologies GmbH
i.V.

Pietro Amatulli
Technical Director


Directive 93/42/EEC, all modifications included



FB 7331	List of safety relevant materials and components	
Revision: A	Identification No.: MMo/14/03/12-01 Identification code/year/month/day number	Page 1 of 4 FB 7331 to CC MSNeXT V07_MMo140312-01.docx

For:

Development project:	MedioStar SE
Project number:	480093A73031
Name of the product (if available):	MedioStar (Brand MedioStar NeXT)
Item number (if available):	1590x, 1593x, 1594x (x = 0 – 9)

Prepared:

Vs.	Change short description	Date	Name	Sign
07	Added to No. 7.3 - new enclosure material (Bayblend T65 + Allruna W55 E203, new manufacturing process)	Mar. 12, 2014	M. Moller	
06	Pt 1.6 + 5.1 PS with 3 rd edition approval	Oct. 10, 2013	M. Moller	Sgd. Moller
05	Added No. 7.2/3 - new enclosure Plastic Handpiece cover 400106103 and Plastic device cover 159506027	June, 20 2013	M. Moller	Sgd. Moller
04	Added new power supply on SOC (MP05-xx, ME05-xx von Cotex) deleted AC/DC Power Module TMT10124 (Traco Power), VCP15US24 (XP-Power)	May 30, 2013	M. Moller	Sgd. Moller
03	Added new enclosure material + change in order number for micro pressure switch	Feb. 21, 2013	M. Moller	Sgd. Moller
02	Added UL labeling for pcb, wiring	May 19, 2011	M. Moller	Sgd. Moller
01	First edition	Feb. 25, 2011	M. Moller	Sgd. Moller

Check and release:	Division	Date	Name	Sign
Checked by	QA (KED)	11.03.14	C. Müller	
Released by	QA	17.03.14	Thieme	

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FB 7331	List of safety relevant materials and components	Asclepiion rcart technology
Revision: A	Identification No.: MMo/14/03/12-01 Identification code/year/month/day number	Page 2 of 4 FB 7331 LoCC MSNeXT V07_MMo140312-01.docx

no	type and item number	manufacturer	operating ratings	nominal ratings	approval	place of assembly
1	Electrical supply					
1.1	Line cord 550704013	Feller Burklin OHG	240 V 8A	250 V 16 A	VDE, CE, UL, CSA, ENEC	outside
1.2	Could apparatus plug 550206080	Schurter / O Heil	240 V 7 A	250 V 16 A	VDE, CE, CSA, UL	inside
1.3	Line fuse/Line switch 550901034	ETA	240 V 7 A	250 V 16 A	VDE, CCC, CSA, UL	inside
1.4	Line Filter B 551201026	ALT/FEG	240 V 7 A	250 V 16 A	Tested by ALT LinefilterB FB7353_Prufu ng Bauteile Komponenten pdf RT/11/02/21-01	inside
1.5	Distribution line 1590 01 017					inside
1.5.1	Terminal Block 788-304 550103055	WAGO	240 V 3A/ 4kV	250 V 3A/ 4kV	KEMA, UL, CSA,VDE	Part of pos. 1.5
1.5.2	Relais RT314024 /88-154 550103056	Tyco/Schrack	240 V 8A/ 5kV	250 V 16A/ 5kV	VDE, UL, CSA	Part of pos. 1.5
1.6	Medical power supply RL0402E-24 550399112	Gunter Dienstleistungen / Ultra Level Tech	In 100- 240Vac ~600 W 47-63Hz Out 24Vdc 350W	In 100- 240Vac ~700 W 47-63Hz Out 22- 30Vdc 400W	EN 60601-1:2006, CAN/CSA-C22 2 No 60601-1-08, ANSI/AAMI ES60601-1:2005 UL E257331 CB 60601-1 (DE 02015163) RoHS	inside
2	Switch-on circuit					
2.1	PCB SOC 7601 06 032					inside
2.1.1	MP05-xx ME05-xx	Cotek	Input 230 VAC, 0.5A Output DC 24V, 0.21A	Input 85- 264 VAC, 0.13A Output DC 24V, 0.21A	UL60601-1 E324816 CUL 60601-1 CB TUV IEC 60601-1	Part of pos. 2.1 PCB SOC
2.1.2	Fuseholder 0031.8201, 0853.0551	Schurter	240 V 6 A	250 V 10 A	IEC60127, VDE, CSA, UL	Part of pos. 2.1 PCB SOC
2.2	Key Switch / Laser Emergency Switch Schaltelement Oktron BZOI 482 550101039 BZOO/1.000.7706 550101040	Schlegel	12 V DC 1 A	250 VAC 6 A	EN 61058 IEC60947-5-1 TUV RU, CSA, UL	inside

FB 7331	List of safety relevant materials and components	
Revision: A	Identification No.: MMo/14/03/12-01 Identification code/year/month/day number	Page 3 of 4 FB 7331 : oCC MSNeX1 V07 MMo140312 01.docx

no	type and item number	manufacturer	operating ratings	nominal ratings	approval	place of assembly
2.3.	Foot switch KF1SD1S MED SK12 660806008	Steute	10 VAC 4 20 mA	max 25 VAC max. 1 A	DIN EN + IEC 60601-1 DIN EN + IEC 60601-2-2 DIN EN + IEC 60601-2-22 DIN EN + IEC 60529 IPX8 UL 60601-1 MPG 93/42/EEC UL 94-V0/-V2 RWTÜV	outside
3.	Control system					
3.1.	ECD4 PCB 760106027	ALT/ACD	+15 V < 2 A	+8-40 V 4 A	UL, Verified by ALT	inside
3.1.1	Relais ALD 105 Relais ALD 112	Panasonic	+5 V < 40 mA 1.5 kV	277VAC/3 A 30VDC/3A 4 kV	UL, CSA, VDE EN 61810-1	Part of pos. 3.1
3.1.2.	DC-DC Converter NMV0512	Murata	1.5 KV	3 kV	Tested by SLG	Part of pos. 3.1
3.1.3	Software 1590(1593)04400-402	ALT	N.A.	N.A.	Tested by ALT	Part of pos. 3.1
4.	Cooling system					
4.1	Micro Pressure switch 580067 Alternativ: 18575: HR30211B3 HP 19732: HR10611B3 HP	Wilspec	12 - 15 Bar	17,7 Bar	UL TUV EN 12263, 97/23/EG (PED) DIN 32733	Part of pos. 4 (inside)
5.	Laser power supply					
5.1	Medical power supply ME 1200-24 550399111	Gunter Dienstleistung n / Codec	In 100- 240Vac 47-63 Hz ~1600W Out 24Vdc/ ~1200W	In 100- 240Vac 47-63 Hz ~650W Out 24Vdc/ ~500W	EN 60601-1:2006 (TA5022871701) CAN/CSA-C22.2 No. 60601-1:08. ANSI/AAMI ES60601-1:2005 UL E324816	inside
6.	Wiring harness					
6.1.	Braid Wire H07V-K 550701117 550701118 550701121 550701122 550701144 550701145	U.I. Lapp	240V	700V	Tested by SLG	inside

FB 7331	List of safety relevant materials and components	
Revision: A	Identification No.: MMo/14/03/12-01 Identification code/year/month/day number	Page 4 of 4 FB 7331 LoCC MSNeXT V07_MMo140312-01.docx

no	type and item number	manufacturer	operating ratings	nominal ratings	approval	place of assembly
6.2	Braid Wire H05V-K 550701116 550701120 550701147	UI Lapp	240V	500V	Tested by SLG	inside
6.3	Protective tubing PVC 550806040	Alpha Wire	N.A.	N.A.	UL, CSA Tested by SLG	inside
Enclosure						
Plastic device cover 159006001						
7.1	Biresin RG53 FR	SiKa Deutschland	Device inside < 50 °C	110°C	RoHS, UL 94, V-0 at 3 mm RU E204929	outside
Plastic device cover 159006023 + 159506027						
7.2	Styrol/Butadien (brandgeschützt – SB V0)	Kern GmbH	Device inside < 50 °C	75°C	UL 94, V-0, RoHS, RU E83071	outside
Plastic Handpiece cover 400106103 + 400106035						
7.3	Biresin VG230	SiKa Deutschland	Device inside < 50 °C	70°C	tested by ALT, RoHS	outside
	Bayblend T65 Allruna W55 E203	Bayer Allod	Device inside < 50 °C	70°C	tested by ALT, RoHS, REACH	outside

User Manual

THE EPI LAB

by Asclepion

(Version 15901)

Page intentionally left blank

1590 99 015_V1 – June 24, 2013

WARNING

This laser device emits high levels of emission in the infrared range that may cause serious damage to the eye if the device is used improperly. Moreover, further risks are possible.

Therefore, this User Manual must be studied carefully before using the device and all warnings and instructions are to be observed during its operation.

It is of particular importance that all persons in the room in which the laser is operated put on the appropriate eye protection during the active laser operation.

Page intentionally left blank

1590 99 015_V1 – June 24, 2013

Table of contents

1	Introduction	1
1.1	General information	1
1.2	Symbols in the User Manual.....	1
1.3	Intended use.....	2
1.4	Copyright	3
2	Safety	4
2.1	General information	4
2.2	Electrical hazards	5
2.3	Biological hazards	5
2.4	Radio interference	5
2.5	Laser	6
2.5.1	Optical hazards caused by laser emission	6
2.5.2	Laser-induced fire hazard	8
3	Technical data	9
3.1	Device models	9
3.2	Specifications	9
3.3	Technical description	11
3.4	Safety units.....	12
4	Start up	15
4.1	Scope of delivery	15
4.2	Installation	16
4.2.1	Device design	16
4.2.2	Labels at the device	17
4.2.3	General information about installation.....	21
4.2.4	Door interlock connection	21
4.2.5	Foot switch	22
4.2.6	Handpiece with integrated skin cooling.....	23
4.2.7	Mains connection.....	24
4.3	Switching on	24

5	Operation of the device	26
5.1	General information on operation	26
5.2	Setup	27
5.3	Setting parameters and starting the treatment	28
5.3.1	Parameters	30
5.3.2	Special function/Tools	30
5.3.3	Activating the laser ready mode and releasing the emission of the laser beam	31
5.3.4	Changing the handpiece	33
5.3.5	Switching off the device	33
6	Application	34
6.1	General instructions	34
6.2	General treatment information	35
6.2.1	Interaction of laser beam and tissue	36
6.2.2	Classification of skin types	38
6.2.3	Overview on zones sensitive to pain	38
6.2.4	Application parameters	39
6.3	Typical indications – overview	41
6.4	Epilation	41
6.4.1	The target: anatomy of hair	41
6.4.2	The phases of hair growth	42
6.4.3	Selection of patients	43
6.4.4	Contraindications	44
6.4.5	Possible side effects	44
6.4.6	Preparing laser epilation	44
6.4.7	Test treatment	45
6.4.8	Treatment	45
6.4.9	Post-op treatment	46
6.4.10	Follow-up treatment procedures	46
6.4.11	Expected Treatment Outcome	47
6.5	Vascular treatment	49
6.5.1	Selection of patients	49
6.5.2	Contraindications	50
6.5.3	Possible side effects	50
6.5.4	Preparing the laser treatment	51
6.5.5	Test treatment	51

1590 99 015 V1 – June 24, 2013

6.5.6	Treatment procedure	52
6.5.7	Post-treatment	52
6.5.8	Follow-up treatments	53
6.5.9	Expected Treatment Outcome	53
6.6	Treatment of acne	54
	Expected Treatment Outcome	55
6.7	Skin rejuvenation	56
	Expected Treatment Outcome	57
7	Accessories	58
8	Cleaning, disinfection and sterilization	59
8.1	General rules	59
8.2	General procedure for cleaning and disinfecting the handpieces	60
8.3	Special procedure for cleaning and disinfecting the handpieces	61
9	Fault messages and fault removal	63
9.1	General information	63
9.2	Fault messages and fault removal	64
10	Customer service	66
11	Regular maintenance, safety checks and calibration	67
11.1	Routine maintenance	67
11.2	Regular safety check	68
12	Disposal	69
A	Appendix: EMC manufacturer's declaration	70

Page intentionally left blank

1590 99 015_V1 – June 24, 2013

1 Introduction

1.1 General information

- This User Manual contains the information required for the intended use of the device.
- It is an integral part of the complete product according to the applicable national, European and international directives and laws on product liability.
- The most important objective is to protect persons against risk situations and the device against damages caused by improper use.
- Regular maintenance according to the instructions given in this manual increases the precision and operability of the device during its entire service life.
- In case of malfunctions or errors that go beyond the scope of intervention of the User contact the Technical Service of Asclepion.

Knowledge of this User Manual is required for the operation of the devices. Therefore, please make yourself familiar with the contents and pay special attention to information concerning the safe operation of the device.

Any deviation from the procedures described herein is prohibited and leads to the expiry of the warranty.

1.2 Symbols in the User Manual

In this manual, the following symbols are used to draw your attention to dangers or to give instructions on operating the device:

DANGER

This symbol indicates a direct danger for the life and health of persons.

Non-compliance with these instructions causes serious danger to health up to life-threatening injuries or even death.

WARNING

This symbol indicates a possible danger for the life and health of persons.


Non-compliance with these instructions can cause serious danger to health up to life-threatening injuries.

CAUTION

This symbol indicates a possible dangerous situation.

Non-compliance with these instructions can cause injuries.

Important information

This symbol indicates important information about the proper use of the device. Non-compliance with these instructions can cause damage to property or malfunctions of the device or in its surroundings.

Tip

This symbol indicates practical tips and particularly useful information that enables the User to use all the features of the device in an optimal manner.

1.3 Intended use

Recommended use

The Epi Lab is a pulsed diode laser (Class 4 laser product) that is used for the removal of unwanted hair, for the treatment of acne, for vascular treatments and for skin rejuvenation. Melanin and blood relatively strongly absorb the range of wavelengths of the radiation emitted by this laser, whereas the surrounding skin absorbs only a little bit of it. In this way, the radiation thermally destroys the target structure while protecting the surrounding skin.

Please refer to chapter 6 *Application* for detailed information about the use of this device.

WARNING

The incorrect using of this device can cause serious danger to health. Therefore, all warnings, instructions and informations of this User Manual are to be observed under all circumstances! The non-compliance is considered as unintended use of the device!
Any modification of this medical device is not allowed !

Contraindications

CAUTION

The intensive laser light in the infrared range can cause thermal damage of skin structures. The contraindications listed below are based on clinical experience in operating laser devices for said purposes of during about 15 years. They do not claim to completeness or unlimited validity.

- Light sensitivity
- Usage of preparations that increase sensitivity to light
- Cellulite therapies in the treatment area (wait a few weeks)
- Cancerous and pre-cancerous lesions in the treatment area
- Haired nevi in the treatment area
- Herpes simplex in the treatment area
- History of Keloids

WARNING

The intensive laser light in the infrared range can cause thermal damage to the eyes. Therefore it is **NOT PERMITTED** to laser the area of the eyes! Occasionally, eye inflammation and a reduction of the visual function were observed although the eyes were protected by eyeshields!

WARNING

The intensive laser light in the infrared range can cause thermal damage to skin structures. Therefore it is **NOT PERMITTED** to use laser treatment for pregnant women or children because clinical studies do not exist for these patients.

CAUTION

Patients with tanned skin – especially fresh tanned - develop more intensive side effects, particularly subsequent pigment changes. They should postpone the laser epilation!
Tattoos in the treatment area cause stronger side effects, too. Moreover, the color will be partly removed. Therefore, tattoos should not be treated.

Contraindications particularly for the treatment of spider leg veins

- Superficial varices with a diameter larger than 3 mm
- Branch and side branch varicose veins
- Ulcer
- Diabetes
- Vasculitis
- Stasis dermatitis
- Podedema
- Neuropathy of the legs
- Corticosteroids treatment in the last three months
- Thrombophlebitis

Restriction of use

WARNING

The laser emission of the device can cause thermal changes of tissue structures.

Treatment may not be performed by anyone other than a practicing physician or under the responsibility of a practicing physician (hereinafter referred to as user for simplification).

Active medical products may only be applied by Users who are adequately qualified to do so, based on their professional training background, expert knowledge and practical experience.

It is expressly pointed out that the use of the device is only permitted for Users that have been instructed to operate it and confirmed the training course by their signature in the Medical Device Logbook.

Among other institutions, Asclepion Laser Technologies also offers courses for laser expertise in medical applications.

Medical Cassification

The Epi Lab meets the requirements of the EC Medical Device Directive 93/42/EEC and of its German equivalent in form of the German Medical Product Act (MPG).

Device type acc. to RL 93/42/EEC II b

UMDNS no.. 16-948 laser, therapeutic

1.4 Copyright

The specifications are subject to change due to further technical development. Contact the Technical Service of your distributor to get latest information.

© Unless expressly authorized, the forwarding and duplication of this document, and the utilization and disclosure of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

2 Safety

2.1 General information

The Epi Lab complies with the requirements of the EC Medical Device Directive (93/42/EEC) and the German Medical Product Act (MPG).

The Epi Lab is a Class IIb medical device as per the directive in Annex IX (i.e. a surgical invasive product that is intended to develop a biological effect – that must not be confused with the laser class!).

Like each device, The Epi Lab may lead to potential risks when operated. The user should be aware of them before starting the system. Among these risks are optical, electrical and biological dangers as well as fire hazards.

Observe the national regulations applicable for the operation of a medical product.

According to the EC Medical Device Directive (93/42/EEC) the users in the European Union is obliged to keep a Medical Device Logbook. A corresponding sample is delivered with the device.

WARNING

The laser device emits high levels of radiation in the infrared range that may cause serious damage to the eye if the device is used improperly.

Absolutely observe the instructions given in this User Manual for equipment and procedures for operating, maintaining, checking or calibrating purposes.

WARNING

The laser radiation of the device can cause, whenever low, the risk of fire. Under no circumstances the device may be used in an explosive atmosphere (classified as AP and APG under IEC 601-1). Take care to preclude the existence of vapors from solvents or flammable liquids (these have been possibly used for cleaning or disinfection) in the room before you start working with the laser.

WARNING

Please note that this product is subject to scheduled technical safety tests that must be performed at annual intervals to ensure the safety for the staff and patients. If the tests are not performed regularly serious health problems can be caused.

The results of such safety tests are to be reported in the Medical Device Logbook.

WARNING

Work for service, repair or modification may not be performed by anyone other than Asclepion Laser Technologies personnel or other bodies duly authorized to handle such work to ensure the safety for the staff and patients.

If unauthorized persons carry out service works, repairs or modifications at the device serious health problems can be caused.

2.2 Electrical hazards

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock! Do not open the device, if you are neither trained nor authorized to do so.

The device is only to be operated with reliably grounded and perfectly installed sockets with earthing contact.

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock!

The device is designed to prevent the penetration of liquids during normal use. If liquid penetrates into the inner system space or the handpiece in exceptional cases, you must turn power off and pull the main power plug. Please contact our Technical Service.

2.3 Biological hazards

WARNING

The contact surface of the handpiece contacts the patient's skin during the treatment. There is the risk of microorganisms being transferred.

The selected handpiece must be deep cleaned and disinfected directly after each application!

2.4 Radio interference

The device complies with the requirements of the EN 60601-1-2 standard. The system is not affected by electromagnetic interference generated by other devices that conform to the same standard. Moreover, the system does not generate electromagnetic radiation beyond the limiting values indicated in EN 60601-1-2. The manufacturer's declaration about the electromagnetic compatibility according to EN 60601-1-2 is enclosed to this manual.

WARNING

Medical devices are subject to specific precaution measures concerning the electromagnetic compatibility (EMC). The device emits low levels of electromagnetic radiation that may have an effect on devices that do not comply with the standard. Unintended settings at these devices could be the result of such an effect.

Please observe the notes of the manufacturer's EMC declaration in the Appendix when installing and operating the device.

WARNING

The device may be influenced by devices that do not comply with the standard indicated above. The result may be unintended settings. Therefore, switch off cellular phones and similar equipment before putting the laser device into operation.

WARNING

The use of accessories and cables not authorized by the manufacturer may result in an increased interference emission or a decreased immunity of the device with the above mentioned consequences.

Use only accessories and cables that have been authorized by the manufacturer for this result.

WARNING

The device must not be stacked / arranged close to or in combination with other devices. If the device is operated close to or with other devices it must be checked if the device is working properly in this arrangement.

2.5 Laser

2.5.1 Optical hazards caused by laser emission

Lasers are classified in accordance with their potential for danger. The most dangerous class is class 4 (the less dangerous class 1).

WARNING

The device is a laser of class 4. That means that the direct and even scattered laser radiation can cause serious damage to the eye if the device is operated improperly.

The protection measures described in the following must be strictly observed.

It is of particular importance that all persons in the room in which the laser is operated put on the appropriate eye protection during the active laser operation.

For fundamental rules on the handling of laser devices you are referred to the international standard EN 60825-1. It is complemented by national regulations providing general protection from dangerous laser radiation. Their purpose is to protect operating personnel and patients during medical applications.

Please follow specific national regulations that may be valid for your location.

It is highly recommended to follow the rules specified below:

1. A laser safety officer should be appointed in writing who has at least the following responsibilities:
 - monitoring of laser operation,
 - supporting the user on issues of safe operation and necessary protective measures,
 - cooperating with expert personnel on occupational safety in fulfilling their tasks, including the information about important issues of laser radiation protection.
2. Laser area

During operation, the area where the maximum permissible radiation level may be exceeded, i.e. the so-called laser area must be delimited and marked with a laser warning sign. The laser area is the room in which the laser is used. **Warning lights** and a triangular shaped yellow laser warning label must be provided at all access points.

Figure 1. Laser warning label (IEC 60825-1)

3. The NOHD (Nominal Ocular Hazard Distance) of this laser system is very great so that the entire laser session room is to be regarded the laser area.
4. Personal eye protection

WARNING

Laser radiation can cause serious damage to the eye.

Everyone present in the laser room during a treatment session must wear laser safety goggles, the patients an eyeshield!

The laser safety goggles and the eyeshield for patients must comply with the specifications defined in the technical data section!

5. Further precautionary measures

- Cover up windows and other openings of the treatment room in order to prevent unwanted emission of laser radiation.
- Restrict access to the treatment room to the patient to be treated and to those persons who will assist and are trained in the handling of this laser device.
- Make sure that personnel trained and assisting in a therapy session know how to shut down the laser in the case of an emergency.
- Remove all metal objects such as watches, rings, bracelets or similar belongings from the working range and refrain from the use of reflecting instruments or other comparable materials where possible.
(Reflecting objects can interrupt the laser beam and guide it to another area than the intended one. Many surfaces, even the ones that seem to be mat, can strongly reflect the emission wavelength of the laser.)
- Point the active laser only onto the area that is intended for treatment.
- Refrain from looking directly into the outlet opening of the handpiece or that of the optical fiber, even if you wear laser safety goggles.
- For phases of non-operation, you should switch to STANDBY mode (in the STANDBY mode accidental release of the laser is prevented).
- Always remove the key from the key holder once power has been turned off and keep the key in a safe place.

6. Responsibility

The user who releases the laser emission is responsible for the laser safety. He/she has to ensure that all preventive measures have been taken **before** activating the laser. It is of particular importance that all persons in the room, in which the laser is operated, including the patients, put on the appropriate laser safety goggles or eyeshield.

2.5.2 Laser-induced fire hazard

WARNING

Because of the laser radiation of the devices there is a risk of fire or explosion if the laser is used in the presence of flammable materials, solutions or gases or in an environment which is enriched with oxygen.

The Users should take the following precautions in order to prevent fires caused by laser radiation:

- Please expose only the proposed body areas to laser emission
- Please use non flammable substances for the preparation of the treatment and for the cleaning and disinfection of the handpieces directly before and during the treatment. If it is not possible to avoid usage of solvents such as alcohol or isopropyl alcohol directly before and during the treatment, the vapour should be removed by aeration before activation of the laser or at least wait for some time.
- Please use non-flammable substances for anesthesia, if at all.
- Avoid the use of oxidizing gases, e.g. nitrogen oxid (N_2O) and oxygen. Be especially careful when using oxygen. Oxygen increases both the intensity and the extent of fire.
- Please keep a minimum of flammable materials in the treatment room only. If it is necessary for the treatment to use flammable materials, please moisten these.
- Keep clothes as far as possible away from the treatment zone.
- Some materials, which are saturated with oxygen (e.g. cotton), could be inflamed at the high temperatures, which develop under the intended use of the laser.
- Have always a small fire extinguisher and water ready in the treatment room

3 Technical data

3.1 Device models

Designation	Description
The Epi Lab (Part no. 15901)	High power diode laser
Handpiece 14 x 10 (Part no. 4060)	Standard handpiece for hair removal
Handpiece 8 x 6 (Part no. 4062)	High-power handpiece for hair removal (optional)
Handpiece 4 x 3 (Part no. 4061)	Handpiece for vascular treatments (optional)

The standard delivery includes the main unit with one standard handpiece.

3.2 Specifications

Specification	The Epi Lab
Device model	Table-top device
Display	8.4" LCD display
Operator guidance	Touchscreen
Cooling	Internal cooling cycle
Door interlock connector	5 V / 10 mA
Laser warning lamp	Potential-free relay contact max. 24 V / 1 A (make contact)
Permissible ambient conditions	Temperature: 15 °C to 30 °C rel. air humidity: max. 70 % (no condensation) height: max. 3000m above sea level
Permissible transport and storage conditions	Upright position (not to lay!) Temperature: 5°C to +40°C rel. air humidity: 10 % to 95 % (no condensation)
Dimensions	365 x 560 x 465 mm (W x D x H) without handles and handpiece holder
Weight	ca. 35 kg
Power requirements	100 -240VAC, 50/60 Hz, max. 1500 VA
Main fuse	overcurrent release 16A, medium slow-blow
Overvoltage category	II (IEC 60664-1)
Max. peak of mains voltage	V MT peak = 2500V peak
Max. nominal mains voltage	V MN r.m.s. = 300V r.m.s.

Classification according to
directive 93/42/EEC

II b

Accessories

see chapter 7

Operating mode

The device has been designed for continuous operation*.

* At high room temperature and high average power it is possible that, after longer treatments, the device goes into the cooling mode for about one minute to guarantee the efficiency of the skin cooling.

Specification

The Epi Lab

Laser model

High-power diode laser array

Wavelength

800 - 950 nm

Laser class

4

Required laser safety goggles
(acc. To DIN EN 207: 2009)

800 - 950 nm D LB5

Pulse duration

max. 400 ms (depending on the mode and handpiece)

Pulse energy

max. 60 J (depending on the mode and handpiece)

Pulse frequency

max. 12 Hz (depending on the mode and handpiece)

Spot sizes (handpieces)

handpiece STD: 14x10 mm²

handpiece HP: 8x6 mm²

handpiece VAS: 4x3 mm²

Tolerance of output power

± 20 %

Laser beam mode

multimode

Laser beam diameter at the
handpiece output

see spot sizes

Energy density

max. 44 J/cm² @ handpiece STD: 14x10 mm²

max. 90 J/cm² @ handpiece HP: 8x6 mm²

max. 210 J/cm² @ handpiece VAS: 4x3 mm²

Beam divergence at the handpiece
output (round angle, 1/e²)

typ. 400 mrad @ handpiece STD

typ. 700 mrad @ handpiece HP

typ. 600 mrad @ handpiece VAS

Nominal ocular hazard distance
(with handpiece)

35 m

Protection class of the part used

B

3.3 Technical description

Power supply unit

The integrated power supply units transform the line voltage in rectified low voltage for the electrical supply of the device.

Control

The output parameters and the safety systems of the device are monitored and controlled by a microcontrol system.

Operating modes

The device operates in both single pulse mode and continuous mode.

The device is designed for continuous operation in the specified ambient conditions. At high room temperature and high average power it is possible that, after longer treatments, the device goes into the cooling mode for a short term to guarantee the efficiency of the skin cooling.

Operation

The device is operated via a touchscreen. The data are output by an LCD with plain-text user guidance. Depending on the operating mode different parameters – such as fluence, pulse duration, pulse rate and others – can be set as described in the Operation chapter.

In addition to the parameters selected, other data such as the laser mode (STANDBY or READY), the handpiece used and, if necessary, error plain-text messages are displayed.

The device contains a database where the user can store treatment reports with the parameters typically used for various types of treatment

Cooling

The internal cooling system is required to prevent the component groups from being overheated. The lost heat is dissipated to the ambience via a water-air heat exchanger. A pump circulates the coolant through the component groups and the heat exchanger. In addition to this, integrated temperature and flow sensors provide safe operation.

The device cooling is ensured by the indoor air. External cooling water is not required

To maintain the cooling effect, the device should only be operated in the specified ambient conditions. The ambient air should be dustless. At the openings provided for cooling the device, a minimum distance should be kept to air-tight objects (e.g. the wall of the room). (See chapter *Installation*).

Laser

The Epi Lab laser is a pulsed high power diode laser (laser class 4).

The emitted wavelength is in the near infrared (NIR) of between 800 nm and 950 nm. The laser beam emission is activated by pressing the foot switch.

Skin Cooling

Next to the laser tip is a metal probe at the handpiece, which is cooled by an integrated Peltier element. Together with the procedure described in chapter *Application* it is guaranteed that the epidermis is cold before the laser pulse impinges. The temperature of the probe is monitored.


3.4 Safety units

The device is equipped with several safety units that have been provided to avoid maloperation and the unintentional activation of the system. All the persons who operate the laser or assist during the treatment should make themselves familiar with these units.

Laser emergency STOP switch

The red Laser emergency STOP switch at the front of the device (see figure below) is for the purpose of switching off the system immediately in an emergency. It should only be activated in emergency situations, i.e. if the emission has to be interrupted immediately.

To turn off the system immediately press this button. To unlock this button after the elimination of the emergency, turn it to the left until it jumps out again.

Note

The Laser emergency STOP switch should not be used to turn the system on and off in its normal condition.

Key switch

The key switch is used to activate the system. Only authorized persons having access to the switch can start the system. The system can only be started via the key switch, if the Laser emergency STOP switch has not been pressed.

WARNING

The laser device emits strong radiation in the infrared range. This radiation can cause serious eye injuries if the device is not properly used.

The device must not be switched on or operated by unauthorized personnel.

Remove the key always after having switched off the system and ensure that it is kept by authorized persons only.


Figure 2: Safety units at the front panel of the device

Power switch

The power switch separates the device from the line voltage. If the device is not in use, the switch should be in the OUT (O) position.

Power switch

Figure 3 Power switch at the rear panel of the device

Foot switch

The foot switch is an electric button that releases the laser emission if the laser is activated and if the device is in the READY mode. It is equipped with two redundant switching elements. Position the foot switch always close to the treatment area.


Figure 4 Foot switch

STANDBY mode

The STANDBY mode prevents the unintentional or inadvertent activation of the laser. If the device is in the STANDBY mode, laser emission cannot be activated. The user can only start the emission, if the READY button has been pressed.

The system switches to STANDBY in the following situations:

- after the first start of the device,
- if the User presses one of the buttons for selecting the laser parameters,
- if the system has been in the READY mode for a longer period of time without the laser released,
- if the User presses the STANDBY key when the system is in the READY mode.

Automatic shut-down

The device is equipped with an automatic shutdown system. If a specific problem arises, the system will automatically change into the following safe mode:

The laser beam shutter is closed, the laser discharge is interrupted and the footswitch is deactivated. The screen shows a fault message identifying the specific error (see chapter 9.).

Remote door interlock contact

The system is equipped with a remote interlock contact that may be connected to the entrances to the laser room. If the remote interlock contact is open (e.g. if the door is opened), the laser emission will be interrupted automatically and the device will change to the STANDBY mode.

Acoustic signal

The laser emission is indicated by an acoustic signal during the period of emission.

Warning labels

The device carries several warning labels (see chapter *Installation - Labels*). These labels must always be clearly readable and replaced immediately when damaged. Contact the Technical Service immediately in such cases.

The manufacturer will only guarantee the safety, reliability and performance of the device, if:

1. the device has been used in accordance with the instructions given in this manual (on the safety measures and the use of the system),
2. the installation, assembly, extensions, changes, repair and maintenance works have been performed by persons being authorized and qualified for them, and
3. the electric system at the place of installation complies with the requirements placed by the EN/IEC and the local regulations.

4 Start up

4.1 Scope of delivery

Note

Check together with the haulier whether all components that are listed below and required for the safe operation of the device are present and not damaged. Incompleteness and/or damages in transit are to be claimed immediately at the supplier.

Asclepion Laser Technologies does not accept any liability unless another agreement has been made expressly.

Please, keep the packing for the future safe consignment of the device.

CAUTION

Refrain from start-up action of any kind if the device or its accessories are found to have suffered mechanical damage. There is risk of physical injury if this advice is disregarded. Contact the Technical Service of your distributor in such cases.

The basic equipment comprises the following device parts:

- basic device
- power cable
- foot switch
- handpiece
- warning lamp/door interlock connector
- keys (two for key switch)
- laser safety goggle, eye protection for patient
- user manual
- laser warning sign

4.2 Installation

4.2.1

Device design


Figure 5: Connections at the front panel of the device


Figure 6: Connections at the rear panel of the device

4.2.2 Labels at the device

WARNING

The laser device emits high levels of radiation in the infrared range that may cause serious damage to the eye if the device is used improperly.

Therefore, warning and information labels are attached at appropriate points of the device. All the staff must know these labels and their meaning to ensure the proper use of the device.

All labels must be kept at their position and must be in good condition. Replace defect labels immediately to ensure the appropriate use of the device.


Figure 7: Labels at The Epi Lab device

Figure 8: Laser emergency STOP switch

This label indicates the Laser emergency STOP switch for switching off the laser radiation promptly.


Figure 9: Beam exit aperture

This label indicates the exit point of the laser beam out of the transmission system.


Figure 10: Interface connections

Warning: It is absolutely important that only the accessories released by Asclepion Laser Technologies are connected here. Master, Service and USB are exclusively designed for authorized personnel.

- 1 Master interface
- 2 USB interface
- 3 Remote interlock contact connector
- 4 Connector for foot switch
- 5 Service interface
- 6 Power inlet socket

Figure 11: Potential equalization connection

This contact is used to connect the device additionally to a potential equalization system.


Figure 12: Type plate The Epi Lab

Type plate indicating type-specific data of the device, such as serial number and date of manufacture.


Figure 13. Laser warning label DIN EN 60825-1

Figure 14. Caution! Consult accompanying documents

Use only coolant liquid "BC4"!

Please note: never mix
BC4 with other liquids!

and

Figure 15. Caution! Observe coolant and upright transportation!

This label informs about the dangers, the maximum values of energy emission and the classification of the laser source.

This label informs the user, that he/she has to read the accompanying documents before the device is used.

This label informs the user, that the coolant has to be used which was specially developed by Asclepion for the device (corrosion protection, antibacterial, frost protection, effective cooling)! This cooling fluid should never be mixed with distilled water.

This label informs the user, that the device has to be kept in a vertical position.

4.2.3 General information about installation

Important information

After the installation and any relocation of the device from a cold to a warm environment with a temperature difference of more than 5°C, allow the device to adjust to the room temperature in unpacked condition before being used at least in the following way:

min. 2 hours at a temperature difference of up to 10°C

min. 4 hours at a temperature difference of up to 15°C

min. 8 hours at a temperature difference of more than 20°C.

Afterwards, switch the device on without a handpiece connected and allow the device to warm up for at least 30 minutes.

Take care that the installation and operating environment of the device meets the following conditions: temperature: 15 to 30°C, relative humidity: max. 70 % and height not more than 3000 m above sea level.

Non-observance of these instructions may result in the destruction of the device.

For an efficient cooling of the device, please keep a minimum distance of at least 20 cm between its lateral and rear sides to the wall.

WARNING

The laser device emits high levels of radiation in the infrared range that may cause serious damage to the eye if the device is used improperly.


Every entrance to the laser area (usually the laser treatment room) must be equipped with warning lamps that are lighting (or flashing) as long as the laser emits radiation to warn persons against entering the room.

4.2.4 Door interlock connection

At the rear panel of the device, a connector is provided for the connection of a door interlock contact to the device (see Figure 16: Connections for remote door interlock contact and laser warning lamp).

As port a round plug has to be used which can be screwed, type (x)V 40 according to IEC 60130-9, e.g. Lumberg WSV 40 (as delivered).

(x): dependent on type, e.g. WS for an angled plug as delivered.

Door interlock contact 5 VDC / 20 mA


Relay contact laser emission (max. 24 V / 0,5 A)

Figure 16: Connections for remote door interlock contact and laser warning lamp

By connecting the two connections *remote door interlock contact* and external door contact you can prevent laser radiation from being emitted when the door of the laser room is opened.

This safety system serves to shut off the surgical beam automatically when the door is being opened thus precluding any risk to the person entering.

The contact connector delivered is factory-fitted with a jumper. To connect a door interlock contact, remove this jumper and connect your door interlock contact to the same terminals. The external contact must be potential-free. If a door-interlock contact is not used, make sure to keep the plug with the factory-fitted jumper connected.

Important information

The contact plug may only be connected with the device being turned off and disconnected from the power supply. Non-observance of this instruction can cause damages of the device.

The contacts for the remote interlock and the foot switch must never be connected with line power because otherwise the system will be seriously damaged. Connect these two connectors only in the way described in this section.

Tip

If you have connected the door interlock contact with the entrance door, please check every day before starting a treatment that the laser cannot be switched to the *READY* mode with this door being open. In case of malfunction please inform the Technical Service.

An external voltage source for a low-voltage warning lamp can be switched via two further contacts of this connector. The internally installed switching relay is designed for a maximum voltage and a maximum current as shown in figure above. Absolutely observe the following safety requirements for the connection of an external warning lamp

WARNING

To avoid the risk of an electric shock with consequential serious health problems as well as fire hazard of the transformer, if you use an external transformer to power the lamp, please note:

A safety transformer complying with IEC/EN61558-2-6 must be selected.

The internal relay of the device is closed, if the device is in the *READY* mode.

The prescribed laser warning lamp may, of course, also be turned on and off via a normal light switch if the laser is operated.

4.2.5 Foot switch

Now, connect the foot switch to the foot switch connector provided at the rear panel of the device (refer to *Figure 6.*). Make sure to insert the plug as deep as it will go. Then, lock it by turning the cap nut clockwise. The foot switch should be positioned right next to the device.

To be able to activate the foot switch later, the cover must be pressed first so that it opens. Thus, the real switch (black) can be accessed and can be activated by the foot to release the laser emission in the *READY* mode.


Figure 17: Foot switch

4.2.6 Handpiece with integrated skin cooling


Figure 18: Handpiece with skin cooling

The laser spot size can be varied by changing the handpiece. The system automatically detects the new handpiece (optional) and considers it accordingly.



Figure 19 Laser tip of handpiece with protective cover

4.2.7 Mains connection

Important information

When selecting the power outlet, consider all technical data of your local network and the device. First, verify that the available line power corresponds with the specified power requirements (line power indicated on type plate). If the device is connected to an inappropriate power supply the device can be damaged.

Please observe the requirements stipulated in your national regulations/standards as amended.

Before the device is connected to the line power supply, the rocker switch of the main switch must be in position $\bigcirc$ (device without potential). First, plug the mains cable into the power inlet socket at the rear panel of the device before connecting the other end with the mains plug. The operation of the laser requires that the device is connected to a single-phase power outlet protected by a separate fuse rated for at least 1500 VA (circuit breaker type C).

4.3 Switching on

Verify that the following requirements have been met before switching on the device:

- The device is in perfect condition (no damages).
- The voltage specified on the type label agrees with the line voltage available at the place of installation.
- The power plug of the device has been plugged into an appropriate power outlet.
- The foot switch plug has been properly connected to the foot switch socket at the rear panel of the device and has been screwed down.
- The plug of the door contact is properly connected and screwed to its terminal at the rear panel of the device.
- Handpiece is connected.
- Warning lamps are switched on at every laser room entrance.

Make sure that all safety precautions have been taken.

Then, proceed as follows:

1. Switch the power switch at the rear panel of the device to I (O).
(Figure 6. *Power switch at the rear panel*)
2. Verify that the Laser emergency STOP is deactivated (not pressed)
3. Then, put the safety key into the key switch and turn it clockwise as far as it will go.

While the device is automatically testing essential and safety-relevant parts, the start screen appears on the display and changes after some time to the treatment screen.

Figure 20. Start screen

5 Operation of the device

WARNING

The laser device emits high levels of radiation in the infrared range.

The use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

This may cause serious damage to the eye.


Absolutely observe the instructions included in this User Manual for equipment and procedures for operating, maintaining, checking or calibrating purposes.

5.1 General information on operation

The treatment parameters are to be set via the touch screen on the display at the front panel. The setting is very straightforward and self-explanatory.

- To select a function press the corresponding fields on the display. By pressing the +/- keys on the right side of the display you can browse through lists consisting of several pages.
- To select parameters press the corresponding parameter field. Fields with a dark background (function fields) can be selected whereas fields with a light background (indication fields) are only used for giving information.
By pressing the +/- keys on the right side of the display you can change the set diameters. If you keep the key pressed for a longer time, the parameters will run through automatically. The new parameters are loaded by pressing OK.
- Successful key pressure will be confirmed by a short beep signal.
- The device is activated by pressing the *READY* key. The active mode is shown in white. If the Device is switched to the *READY* mode, the yellow warning sign will be additionally shown on the display.

Overview on important functions of the device

Function/ parameter	Display	Remarks
Setup		General settings for the device can be made via the display.
Tools		Saved special functions for the laser can be called.
Setting keys		Scrolling through lists and changing settings
Back		Returning to the last screen without loading the changed settings
OK		Loading the changed setting (otherwise the old parameters are activated again).
Activation of footswitch		The active mode is shown in red (here shown STANDBY, the footswitch cannot activate emission). If READY is activated (red) additionally the yellow laser warning sign appears at the right upper corner of the display.
Switching off		The device can be switched off without turning the key switch.

5.2 Setup

You can change the general settings of your display such as language and sound volume. Call the Setup menu by pressing the Setup key in the main menu. The Setup screen appears.


Figure 21 Setup

You can execute the following functions by pressing the corresponding keys:

Function/ parameter	Display	Remarks
Information		Indication of software version
Language		Selection of different languages
Sound		Setting the sound volume
Time lock manager		Input of the code to extend the use time (the code is supplied by your distributor)

5.3 Setting parameters and starting the treatment

After start screen the treatment screen appears after some time.


Figure 22: Treatment screen (example)

On the display you can first select the different modes – Basic, Professional, Smooth Pulse and Vascular – if they are released for the connected handpiece.

Then select the parameter you want to change by pressing the corresponding field. The selected parameter will be visually marked (darker). Use the setting keys +/- to make changes.

The parameters that can be set depend on the mode.

Basic mode:

energy density (Fluence) and treatment speed (Speed) can be set

Note: In this mode, the shortest possible pulse duration (Time) for the set fluence will be used, always with double pulse.

Professional mode:

Fluence, Speed, Time as well as Single or Double pulse can be set

Note: For the double pulse Fluence and Time always refer to the sum of the two parts of the double pulse.

Vascular mode:

Fluence, Speed and Time can be set

Note: always single pulse

Smooth pulse mode:

Fluence and Speed can be set

Note: always single pulse, always shortest possible time for set fluence
low fluence, high frequency
contrary to the other modes, several passes are always required
(see Application Epilation)

Press OK for saving the changed settings (otherwise the old parameters are activated again) !

Tip

Please note that

- only the useful modes can be activated for a specific handpiece. That means that not all modes can be set for each handpiece.
- most of the parameters technically depend on each other. That means that, for example, the highest fluence is not simultaneously possible with the shortest pulse duration.

The spot size can be changed by changing the handpieces and is automatically recognized by the device and displayed.

In addition to the parameters selected, the device status [*STANDBY* or *READY*] and, if applicable, plain-text fault messages are displayed.

Furthermore, the pulse energy, the currently emitted pulse number and the total pulse number are

indicated additionally (in Tools

The laser is released by pressing the foot switch.

5.3.1 Parameters

Function/ parameter	Display	Remarks
Energy density	Fluence	Adjustment of fluence
Treatment speed	Speed	Adjustment of treatment speed (pulse frequency)
Spot size	Spot	Information about the spot size of the handpiece
Pulse duration	Time	Adjustment of the pulse duration <i>Professional</i> : for double pulse sum of the two pulse parts
Operating mode	Mode	Operating mode
Pulse	Pulse	<i>Professional</i> : switching between single pulse/double pulse <i>Burst</i> : number of pulses

5.3.2 Special function/Tools

Press the *Tools* key to open the *Tools* screen of the application).


Figure 23 Tools screen


Function/ parameter	Display	Remarks
Pulse counter	Actual	Current pulse counter
Pulse counter	Total	Total pulse counter
Pulse counter	Remain	Remain pulses till end of life time
Temperature skin cooling	Skin Cooler	cold / medium / warm (recommendation: cold)
Cooling water pump	Pump OFF	Is switched on automatically after changing the handpiece
Energy indication	Display energy	Indication of the energy on the display

Important information

Before the real counter reaches end of life time (and "Remain" reaches 0) there are several reminders on the display that the attached handpiece has to be send in for revision to the service department of Your distributor in a foreseeable time.

This way a long term planning is possible. During revision it is decided, if the handpiece can be used further for a specific number of pulses or if the handpiece has to be replaced.

Important information

Before changing the handpiece the cooling water pump must always be switched off by switching off the device (recommended) or pressing the *Pump off* key in the Tools menu. Thus, leakage of the cooling agent is prevented. When the handpiece has been changed the pump is switched on automatically.

5.3.3 Activating the laser ready mode and releasing the emission of the laser beam

See also chapter *Application*.

After the switching on of the laser device it is in the *STANDBY* mode.

Verify again that

- the handpiece has been connected correctly,
- all laser safety precautions have been taken.
- the door interlock contact functions properly, if you have installed it: Open the door, it must not be possible to switch the laser to *READY* (an error message is displayed and must be confirmed).

Set the desired parameters according to the explanations given in the chapters before.

WARNING

The direct and scattered laser radiation can cause serious eye injuries.

During the active laser operation all persons present in the laser room have to wear protective eye goggles or patient eye shields appropriate for this laser as specified in chapter *Technical Data*. Other types of laser protective eye-wear do possibly not provide the degree of protection needed!

Before pressing the *READY* key the attending user must check this.

In the *READY* mode you can also activate a laser pulse by pressing the foot switch unintentionally! Therefore, the laser device must be switched to *STANDBY* in the treatment pauses and at the end of the treatment.

The laser beam may only be released if the handpiece is precisely directed onto the skin area to be treated.

In case of dangerous malfunctions the Laser emergency STOP must be activated immediately.

By pressing the *READY* key you change to the *READY* mode; the yellow laser warning label appears on the bottom right of the display. First, the system performs a safety test of the hardware that takes about 2 sec (see Figure below).

Moreover the system checks, whether the temperature of the cooling probe of the handpiece is low enough. If the cooling temperature after the hardware start has not been reached yet, the Coldown handpiece screen is displayed. Its progress bar reflects the cooling-down of the cooling probe. This screen will be faded out, if the cooling probe is cold.

(The same applies for activating the device after 15 minutes of inactivity.)


Figure 24 After activation of READY

If the footswitch is pressed, the laser radiation is released from the laser tip of the handpiece. The release of the foot switch instantly stops the emission of laser radiation.

Emission of laser radiation is indicated acoustically by a warning beep and optically by the blinking laser emission indicator on the display.

When you change any parameter in the *READY* mode, the device will automatically switch to *STANDBY*, i.e., you have to press the *READY* key once more to return to the *READY* mode.

During the treatment it is normal that the temperature of the cooling probe of the handpiece heats up a little bit, particularly for larger areas. If the upper limit of the cooling temperature is reached, the emission of the laser radiation will be interrupted and an informational screen will be faded in. The treatment can be continued as soon as the cooling tip is cold enough again.

The Epi Lab is equipped with a user-friendly Laser emergency STOP that is located at the front panel of the device. In an emergency, the laser is immediately deactivated by strongly pressing this key and thus any laser danger for the user and patient is removed.

If the emergency problem has been solved, restore the functionality of the device by slightly turning the Laser STOP switch manually anticlockwise to its normal position.

CAUTION

If during lasing the symbols $\triangle$ or ∇ , appear in the Fluence field, the device has detected a deviation of the output values higher or lower than 20 %.

Higher values can cause undesired side effects, whereas lower values can reduce the effect of the treatment.


If this deviation continues inform the Technical Service.

If you have finished the treatment, switch the laser to the safe *STANDBY* mode by pressing the *STANDBY* key. If you intend to interrupt the treatment for some time, switch the laser to the *STANDBY* mode, too. In this mode, it is impossible to release laser emission unintentionally.

After a longer period of inactivity in the *READY* mode, the device automatically switches to the *STANDBY* mode.

5.3.4 Changing the handpiece

The handpiece shall be principally changed with the device being switched off. If this is not possible, the cooling pump must be deactivated imperatively before the handpiece is removed from the device. For doing this, press the Pump OFF key in the Tools menu to switch off the pump and the LED at the handpiece extinguishes. Only now the handpiece plug can be separated from the device without any risk.

Important information

Before changing the handpiece, the cooling water pump of the device must be switched off by pressing the Pump OFF key in the Tools menu. When the handpiece has been changed, the system switches on the pump automatically.

5.3.5 Switching off the device

Exit the *READY* mode and switch to *STANDBY* mode by pressing the *STANDBY* key.

To switch off the laser device, turn the key switch anticlockwise as far as it will go. Only then, you should switch off the power switch at the rear panel of the device. By switching the power switch to '0' the device is completely disconnected from the line voltage.

WARNING

The direct and scattered laser radiation can cause serious damage to the eye.

Never leave the device unattended while it is switched on. Shut down the device completely by switching the power switch.

After you have turned off the device, make sure you have removed the key from the key switch. Keep the key in a safe place to prevent the unauthorized use of the device.

6 Application

6.1 General instructions

This chapter contains detailed information on the operation and application of the laser device. Of course, the instructions given can never replace the expert knowledge for the medical application of lasers, the careful study of relevant literature, the personal experience gained under supervision and the critical assessment of the current situation.

But these instructions shall be a useful support for beginners and those not working regularly with the laser.

Therefore, we recommend you to study current medical literature and to contact physicians in private practices using such type of equipment in order to familiarize with the methods of laser treatment before you start to treat patients yourself.

We would like to assist you in getting contact with other users. The responsible partner of Your distributor will readily inform you accordingly in detail.

Demands placed on the user

WARNING


The laser radiation of the device can cause thermal changes in the tissue structures.

Treatment may not be performed by anyone other than a practicing physician, or under the responsibility of a practicing physician (referred to as user for simplification hereafter).

Active medical products may only be applied by users who are adequately qualified to do so, based on their professional training background, expert knowledge and practical experience.

Make sure you have understood the principles of laser-tissue interaction, the correlation between the individual application parameters, the application techniques as well as the basics of laser safety.

If you have the slightest doubt, consult one or more colleagues with practical experience and/or application engineers of Your distributor before you start a laser treatment.

Users who want to use this device must have been sufficiently informed about the safe and efficient operation of the laser equipment described in this manual before starting to operate it. Persons employed with or authorized by Asclepion Laser Technologies give instructions that are based on this manual. This is to be reported in the Medical Device Logbook. The Users must get instructions about the following aspects:

- operation and intended use of the device, practical exercises included,
- operating principle and functionality of the device as well as energies used,
- adjustments of all operating elements,
- indications for the use of the device,
- contraindications and side effects,
- warning labels in all operating modes,
- performance of functional tests.

Further training requirements vary in the different countries. The user takes responsibility for the compliance of the training with all requirements of the applicable regional laws and regulations.

6.2 General treatment information

CAUTION

The application of laser radiation on the skin can lead to thermal damages of tissue structures – intentionally according to the indication or unintentionally. The physician is responsible for the precise determination of the indication, the treatment zone and the treatment parameters as well as for the correct positioning of the handpiece. Also during the treatment the temperature of the cooling probe of the handpiece has to be checked by hand as well as the adjusted Fluence at the display.

CAUTION

Due to the effect of the laser beam side effects, such as inflammation and slight oedemas, may be observed in the treated area during the treatment. Blisters with subsequent incrustation and temporary pigmentary changes are seldom. Persistent pigmentary changes, scars and temporary urticarial changes are extremely seldom if the device has been operated as intended and there is the compliance of the patient. The degree of these side effects varies from patient to patient; normally it is very low. Therefore, a test treatment must be performed in any case in a small area as described below to minimize the side effects and to optimize the effect of the treatment simultaneously. It is important to bear in mind that the side effects can be stronger for the treatment of larger areas than for the treatment of a smaller test area.

Besides, the patients get a better impression of the treatment procedure and can better decide for or against a treatment before giving their consent.

Tip

In a pre-operation discussion, the patient should be informed about all risks that may be caused by the laser treatment and the pre- and post-operative care. Moreover, the patient should get all information about the mechanism of laser treatment, alternative procedures, the time required for the treatment, the prospects of successful treatment and about the behavior after the laser treatment. You are advised to report this conversation and to make the patient sign a special consent form for this briefing. It is also recommended to take photos before and after the treatment.

WARNING

The contact surface of the handpiece contacts the patient's skin during the treatment. There is the risk of microorganisms being transferred. The selected handpiece must be disinfected before each application and cleaned and disinfected after each application session as described in chapter *Cleaning and disinfection*.

CAUTION

Any absorbing contamination on the laser tip of the handpiece generates local overheating and can cause increased side effects and defects at the handpiece! Make sure that the laser tip of the handpiece remains also clean during the treatment!

6.2.1 Interaction of laser beam and tissue

The goal especially for dermatological and aesthetic treatments is to remove the lesions to be treated or the unwanted lesions without any side effects. The selective photothermolysis method introduced by Anderson as a treatment concept in 1983 ensures an effective damage of the target and the simultaneous protection of the skin surrounding it.

Therefore, the following aspects are to be observed for a successful treatment:

Absorption of the laser beam – wave length

The absorption of the laser beam must be higher in the target (e.g. the melanin in the hair or the blood in the vessel) than the one in the surrounding skin mainly consisting of water. The higher the absorption in the target the greater the heat generated; the lower the absorption in the surrounding skin area, the lower the unwanted heating of the skin.

The wave length range of The Epi Lab device in the near infrared (IR) between 800nm and 950nm has a relatively high melanin and blood absorption. The wavelengths of the different handpieces are optimal adjusted within this range to the individually planned indication.

However, you have also to consider the absorption of melanin in the epidermis of dark skin that is generally not desired (apart from the treatment of benign pigmented lesions) and can produce side effects like stronger bad sensation, hyper- or hypo-pigmentation. Therefore, the Rubin and the Alexandrit lasers having wave lengths of 694nm and 755nm, respectively, are less suited for treating dark skin types due to their relatively high melanin absorption. The wave length range of The Epi Lab is considerably better suited.


Figure 25: Absorption coefficients depending on the wavelength of the laser beam

Spot size

The spot size must be so large that the target located some millimeters beneath the surface of the skin (e.g. a hair root or vessel) is reached with a sufficiently high energy density (FLUENCE). This requirement is met by The Epi Lab laser device.

It must be principally ensured that a larger spot heats a larger volume in the skin. Therefore, the required fluence is lower for larger spots.

Pulse duration

The pulse duration of the laser should be shorter or at least not significantly longer than the thermal relaxation time of the target (e.g. hair or blood vessel). The Thermal Relaxation Time is the time in which the target gives off 50% of the heat to the environment.

Thin structures (hair, vessels) give off the heat much more rapidly than thick structures. To achieve the desired temperature in thin targets it is necessary to deliver the required energy as quick as possible. If the pulse is too long, the heat loss into the surrounding skin must be compensated by a higher fluence, which causes by tendency more side effects.


Figure 26. Thermal Relaxation Time and pulse duration of the laser beam

Target	Diameter	Therm. Relaxation time
Ectatic vessels	10 μm – 50 μm	50 μsec - 1 msec
	100 μm	5 msec
	300 μm	45 msec
	500 μm	125 msec
	1mm	(500 msec)
Hair bulge	ca. 300 μm	ca. 45 msec
Melanosomen	ca. 0.5 – 1 μm	ca. 0.05 μsec

Figure 27: Thermal relaxation time for different target diameters (examples)

However, for an effective treatment the required fluence and spot size are defined first. Only afterwards the shortest possible pulse is set for thin structures and light skin.

For dark skin it is important to consider that the melanin in the epidermis shall not be heated, i.e., a long pulse duration should be selected – no problem for thick structures.

A compromise must be made when treating thin hair/vessels of patients with dark skin. The pulse to be selected is longer than for light skin and thin target but shorter than for dark skin and thicker structures. This decision is acceptable because the development of heat of thin targets is considerably lower than of thick targets. Moreover, The Epi Lab ensures an effective cooling of the epidermis (see below).

The epidermis can be protected additionally by the pulse interruption (double pulse). In BASIC mode the double pulse is always emitted, in the PROFESSIONAL mode it can be activated on the display.

Cooling the epidermis

To avoid unwanted side-effects, the effective cooling of the epidermis is necessary (particularly important for dark skin types because of their melanin absorption).

The skin cooling system of The Epi Lab is integrated in the handpiece. The contact of the epidermis with the cold metal probe is extremely effective due to the high thermal conductivity (8 times better than a sapphire-based cooling mechanism). Moreover, the cooling temperature of the skin cooling system is automatically controlled at any point of time.

6.2.2 Classification of skin types

Skin type		
I (‘Irish’ skin type)	extremely fair skin, blond or red hair, freckles, blue or green eyes	always sunburn no tanning
II (frequent in Central Europe)	fair skin, blond hair, freckles blue or green eyes	mostly sunburn slight tanning
III	light brown or brown hair, brown eyes	sometimes sunburn good tanning
IV (Mediterranean type)	light brown skin, dark or black hair, brown eyes	seldom sunburn fast tanning
V (some Indians, South Americans)	dark skin, black hair, dark eyes	seldom sunburn very fast tanning
VI (some African, Indians)	very dark or black skin, black hair, black eyes	never sunburn, very fast tanning

Figure 28 Classification of skin types according to Fitzpatrick

6.2.3 Overview on zones sensitive to pain

Face


Figure 29: Face zones sensitive to pain

Body


Figure 30: Body zone sensitive to pain

6.2.4 Application parameters

The energy density (FLUENCE) is the crucial parameter for the treatment with pulsed lasers like The Epi Lab. It is to set directly on the display. The higher the fluence the more efficient the treatment – but stronger side effects are possible. Therefore, test treatments are to be performed if the fluence is increased.

The mathematical relationships are as follows:

$$\begin{aligned} \text{Energy density [Joule/cm}^2\text{]} &= \text{energy [J] / spot-area [cm}^2\text{]} (\text{energy density} = \text{Fluence}) \\ \text{Power [Watt]} &= \text{energy [J] / pulse duration [sec]} \end{aligned}$$

The spot size used (depending on the handpiece) is automatically considered by the device. As described above, smaller spots require a larger fluence.

The pulse duration can be set in the *Professional*, *Vascular* and *Burst* modes according to the explanation given above. shorter (with sufficient fluence) for thin hair and vessels as well as for pigmented lesions and longer for protecting dark skin. Pauses between the pulse-parts (double pulse, *Burst*) are also a more gentle treatment for the melanin in the epidermis.

The *Basic* and *Professional* modes are intended to be used for epilation. The *Vascular* mode is provided for vascular treatments and benign pigmented lesions (with the small spot of the VAS handpiece). The *Burst* mode can be used for vascular treatments, too.

The *Smooth Pulse* mode can also be used for the epilation of larger areas. Contrary to all the other treatments, in this mode the handpiece with the set pulse sequence (speed) and low fluence is moved rapidly over the area to be treated. However, several passes are required to be able to see a perifollicular erythema.

To ensure sufficient cooling of the epidermis, the handpiece is always drawn towards the cold metal probe in all treatments (apart from pigmented lesions) so that the laser pulse always reaches the pre-cooled skin. The only exception is the treatment of pigmented lesions because the melanin in the epidermis shall be destroyed.

The set pulse sequence (speed) in Hz (pulse per second) determines the treatment speed and does not influence the effect of the treatment. In sensitive areas a reduction can be less painful.

6.3 Typical indications – overview

- epilation
- vascular treatments
- benign pigmented lesions

The following sections contain treatment instructions for the individual indications.

6.4 Epilation

The removal of unwanted hair growth is undoubtedly one of the oldest cosmetic treatments performed by men. Over the past years, it has experienced greatly increased popularity. The methods range from regularly repeated shavings to a more lasting, though very painful wax epilation, to extremely time-consuming hair removal by means of electrolysis needles, which on the other hand can yield a rather lasting success.

Each of these methods, and, naturally, also laser epilation which is comparatively more pleasant, aims to remove hair with a maximum lasting effect.

6.4.1 The target: anatomy of hair

As everybody knows from one's own experience, lasting removal of hair cannot be achieved by simply plucking out unwanted hair. It will grow again. This means that the structures responsible for hair growth have to be destroyed, too.

Today researchers are virtually unanimous that the so called hair bulge with the germinative cells around it (matrix cells, melanin) and the papilla contribute to the hair growth. Therefore, they must be destroyed or at least impaired in their action for this reason.


Figure 31. Anatomy of the hair follicle

Depending on the body part, the hair's roots may be located down to 5 mm beneath the skin surface. Therefore, the laser spot should be as large as possible to ensure a depth effectiveness of laser radiation.

A clinical study performed with The Epi Lab device showed that a spot size of 1.1 cm² is sufficient to treat hair on the legs effectively. A further enlargement of the spot did not achieve better results. The smaller spot size of 0.5 cm² led to a lower effectiveness of ca. 10%.

A point to bear in mind is that, for example at the upper lip, the hair root is not located very deep in the skin and therefore a smaller spot size can be absolutely effective and useful. The use of smaller spots is recommended particularly for finer hair because shorter pulses can be set with the smaller spot.

6.4.2 The phases of hair growth

In the embryonic state hair is already growing – but it is very fine up to infancy. This lanugo does not have a papilla yet, grows from the lower epidermis and successively migrates into the dermis.

Later on, terminal hair with papilla forms under the influence of hormones. The number of hair is dependent on genetics. From this point on, the human hair is subject to repeating growth cycles which may extend from some months to some years (depending on the body region).

The active growth phase is termed anagen, whereas telogen is the phase of resting. To achieve efficient laser treatment, the transition phase to re-growth appears to be particularly suitable.


Figure 32: Hair growth cycle

In humans, the phases of different hair are not synchronized, i.e., different hair is in different phases. In a certain area only a specific percentage of hair reacts on the laser treatment. Therefore, a series of laser sessions is always necessary.

The shorter the telogen phase compared to the anagen phase the more rapid is the growth of the hair in a specific area. That's why it is necessary to shave the beard area more often than other parts of the body. Consequently, the interval between the subsequent laser treatments must be shorter, too.

	Depth of terminal anagen follicle	Telogen duration	Anagen duration	Follicle no. per cm ²	Total follicle number 1 Million	Growth per day
HEAD						
Upper lip	1-2.5 mm	6 weeks	4 months	500		
Chin - beard	2-4 mm	10 weeks	1 year	500		0.38 mm
Cheeks	2-4 mm			880		0.32 mm
Eye brows	2-2.5 mm	3 months	4-8 weeks			0.16 mm
Ear		3 months	4-8 weeks			
Head hair	3-5 mm	3-4 months	2-6 years	350		0.35 mm
BODY						
Armpits	3.5-4.5 mm	3 months	4 months	65		0.3 mm
Legs	2.5-4 mm	5 months	4 months	60	370 000	0.21 mm
Pubic region	3.5-5 mm	3 months	4 months	70		
Arms		4 months	3 months	80	220 000	0.3 mm
Breast	3-4.5 mm			65		0.35 mm
Torso	2-4.5 mm			70	425 000	0.3 mm

Figure 33: Hair growth (extract from Richards-Meharg Table)

6.4.3 Selection of patients

Application of laser has been shown to achieve an extended delay in hair growth both in women and men.

Laser epilation of dark hair is particularly effective because dark hair has higher melanin content. Light hair typically requires higher energy densities that are tolerated because of the normally lighter skin. However, the success of treatment cannot be predicted in each particular case. Surprisingly good results have even been observed with fair hair. On the other hand, there are negative exceptions of particularly resistant dark hair.

Very fair hair, reddish hair (contain phaeomelanin) and grey hair react very poorly or not at all to a laser treatment. The same applies for very thin hair.

Before treating women with a very strong hair growth, the hormonal status should be controlled.

The hormonal status of transsexual men (man to woman) should also be controlled before starting a treatment.

A complete hair removal cannot be achieved, normally fine light hair cannot be removed. It is important to inform the patients about these facts!

6.4.4 Contraindications

For contraindications please see chapter *Intended Use* (besides for leg veins).

CAUTION

Patients with tanned skin – especially fresh tanned – can develop more intensive side effects, particularly subsequent pigment changes. They should postpone the laser epilation!

Tattoos in the treatment area cause stronger side effects, too. Moreover, the color will be partly removed. Therefore, tattoos should not be treated.

A delay of 4-6 weeks is also recommended when other epilating procedures, such as waxing or plucking, have been used because the target chromophore, the 'target' is not available any longer and therefore the laser treatment does not have an effect.

6.4.5 Possible side effects

All side effects are less significant for light skin and light hair and become more intensive for darker, thicker and dense hair and for darker skin and, of course, in case of too high fluence for the individual case (test treatment!).

Therefore, the start values for test treatments are to be selected lower the darker the skin (and normally also the hair color).

During the treatment the patients have a pricking feeling connected with heat.

After the treatment a slight, perifollicular erythema can be seen. If the hair is dark and thick the formation of a perifollicular oedema is also possible. These effects are normal and last for a period of between few hours and 2 days.

If the device is used as intended and there is the compliance of the patient blisters, crusts as well as pigment changes seldom develop (light protection!). Normally, blisters and crust heal without any consequences within about one week. Pigment changes can be observed for some months.

Persistent pigment changes and scars are very seldom.

In rare cases dark hair can decolor and then it cannot be removed by the laser/light due to the missing absorption. It is also very seldom that a paradox hair growth, i.e. stimulation of hair growth at the border of the treated area, is observed.

6.4.6 Preparing laser epilation

As a compulsory pre-condition, the area to be treated must be absolutely shaved and cleaned. Remove also deodorants or other cosmetic products.

Before treating the upper lip a moistened piece of gauze or the like should be placed between upper lip and teeth, in order to intercept the small part of the laser radiation, which can pass through the upper lip, and to avoid the heating-up of the teeth (can be possibly painful).

To document the treatment success photos should be taken.

Normally, anesthesia is not required because the laser epilation with The Epi Lab causes relatively little pain thanks to the precooling with the cold metal probe of the handpiece. For sensitive body areas or very extensive treatment areas, a topical anesthetic, for example EMLA, can be applied depending on the pain tolerance and the application technique (skin cooling!).

A colorless ultrasound gel optimizes the laser beam coupling into the skin to reduce the reflectance and scattering of the laser beam. Thus, a higher efficiency is reached.

The application technique used for the test treatment is also to be used for the later real treatment. A thin layer of clear ultrasonic gel is to be applied on the area to be treated. This also reduces the levels of emerging smoke/smell.

6.4.7 Test treatment

If the patient decides in favor of a laser epilation after getting the information from the physician, a test treatment should always be performed in the area that is to be treated later. The patients' individual skin reaction and pain sensitivity can vary.

The darker the patient's skin the lower the start fluence and the longer the pulse duration to be selected in *Professional* mode (including double pulse). The latter is not a problem for thicker hair. For thin hair the pulse must not be as long so that the efficiency of the treatment can be ensured. Due to the low absorption 'mass' the side effects are normally lower for thin hair.

The application technique used for the test treatment is also to be used for the later real treatment. For doing this, at least 4-6 test spots of the test FLUENCE in *Basic or Professional* mode are to set as a quadrate, because a single spot can cause fewer reactions (better cooling by the surrounding tissue). With *Smooth Pulse* mode an area of 10 x 10cm² is treated.

WARNING

Laser radiation can cause serious eye injuries.

Everyone present in the laser room during a treatment session must wear laser safety goggles including eye protection for the patient according to the instructions given in chapter 3 *Technical Data*!

The set fluence as well as the temperature of the skin cooling probe must be controlled before and during the treatment!

The Epi Lab is switched to the *READY* mode with the desired FLUENCE.

At the start of the operation, the cold probe of the handpiece is shortly placed on the first point to be treated. Then, the handpiece is moved by the spot size so that the laser tip comes to the pre-cooled area. Afterwards, the laser emission is started by pressing the foot switch and the handpiece is further moved by the spot in the cycle of the set SPEED in direction to the skin cooling probe.

Practically, the handpiece should be kept in contact with the skin and should not always be put on anew. Due to the ultrasound gel the point can easily be seen at which the spot must be placed for the next line. A slight overlapping is necessary so that no hair line is left out because the cooling probe, which marks the track in the ultrasound gel, is a little bit larger than the active laser spot.

Only if a skin reaction cannot be seen after 10 minutes the energy density (FLUENCE) can be increased gradually (test on a new area).

In the *SmoothPulse* mode – contrary to the procedure described above – an area of 10x10cm² is pre-cooled and then the handpiece is moved rapidly (ca. 10 cm per second) over an area of about 10 x 10 cm² in the cycle of the set speed. In doing so the handpiece has to be held in a way, that the laser spot moves over the pre-cooled line and the cold probe pre-cools the next line simultaneously. The handpiece can be moved then in a circle or serpentine pattern back and forth. Several passes must be performed.

Make sure that the cold probe of the skin cooling is kept in contact with the skin during the treatment to guarantee the pre-cooling of the epidermis.

It is advantageous for epilation to apply a slight pressure by the handpiece onto the skin. Thus, the distance to the hair root is reduced and the effectiveness of the treatment is increased.

If the area is completely lasered, the device must be switched to the *STANDBY* mode.

6.4.8 Treatment

Once the treatment area has been prepared, the area can be marked by using a white pencil (e.g. Kajal). This ensures that the treatment area will be lasered evenly.

The FLUENCE that was well tolerated in the test treatment is to be set in the real operation. Then, the area is scanned as described above and the handpiece is always moved towards the skin cooling.

In the *SmoothPulse* mode the number of passes must be increased step by step.

The simplest method to keep the tip also clean during the treatment is the use of a slightly moistened cloth after having lased one line. But if a crack can be heard when treating dark, dense hair, one hair rest is mostly on the tip. The tip must be cleaned immediately. If the laser treatment is continued nevertheless, the hair will burn in and can possibly not be removed by wiping it off.

CAUTION

Dirt settling on the contact surface of the laser tip may burn in when exposed to laser radiation and cause side effects such as reddening, possibly blistering, incrustation, seldom scars!

Ensure that all accessible optical surfaces are not damaged and remain clean during the treatment.

If large areas need to be lased, the skin surface should preferably be treated in longitudinal direction (apart from the *SmoothPulse* mode) in order to keep effects of physical pain building up at a minimum level – particularly for dark, thick and dense hair (use EMLA if necessary). Once the area is completed, the laser must be switched to STANDBY mode.

6.4.9 Post-op treatment

Immediately after completion of lasing, the skin should be cooled over an extended period, particularly for dark hair and dark skin. Continuous cooling can considerably reduce the side effects! This can easily be done with CoolPacks which the patient can apply as appropriate.

An anti-inflammatory ointment may be applied following sufficient cooling of the skin. In the event of overtreatment, a topical antibiotic is indicated and the monitoring of healing is absolutely necessary.

CAUTION

If the treated area is exposed to sun and solar radiation after the treatment, stronger side effects can develop, pigmentary changes in particular.

Therefore, such a radiation should be avoided for a period of 6 weeks or a sun protection cream with a protection factor not below 20 should be used at least.

The lased area should not be traumatized either for 2 to 3 days.

Exfoliation can be supported by using a lotion. Patients should be informed that up to two days after treatment the openings of hair follicles may secrete tissue liquid. Possibly resulting crusting should be removed by normal washing and application of skin lotions.

6.4.10 Follow-up treatment procedures

The decision about further treatment procedures should consider the patient's individual skin reaction and his/her personal interests, and the energy densities should be selected to be optimally suited for the given purpose.

The patients are to be informed that in some cases – particularly for thick hair – it takes two weeks that the hair is pushed out.

If the fluence of the last session was well tolerated, it can be increased gradually in the follow-up treatments. This increase is also recommended because the hair becomes less and less dense and thinner. Consequently, a higher fluence is required and normally well tolerated.

The number of required treatments depends on several factors:

- genetic disposition
- hormonal status
- the body region to be treated
(different number of hair in the anagen phase, see *Extract from Richards-Mehharg table*)
- hair and skin type
- result desired by the patient (often, remaining fine hair are accepted)
- interval between the treatments.

To achieve a satisfactory result, four to seven treatment sessions may be necessary on average. Especially for the beard area and for people with very strong hair growth the number of sessions required is even higher.

The same goes for men with back or breast hair that grew only shortly ago. Then, it often seems that the laser treatment does not have any effect but with no treatment the hair growth would have been stronger. In this case it is helpful to ask the patient about the father's or uncle's hair growth at this age

The intervals between the sessions should not be too short (4-6 weeks; in case no hair regenerated the time interval can be longer) in order to allow as much hair as possible to change from the rest to the anagen phase.

In general, the interval for treating the beard area should be shorter (ca. 4 weeks) and for treating the other areas longer (ca. 6 weeks) as explained in chapter 6.4.2.

Normally, the interval is selected longer after several sessions because less hair is growing anew.

However, if a strong growth should have begun after the laser treatment, it is advisable to repeat the treatment in short time intervals because then very much hair is in the anagen phase and shows a good reaction to the treatment

For patients with normal hair growth the achieved result can be stable for several years. For patients with genetically caused strong hair growth, typical for people of southern countries, it is to be expected that a follow-up treatment is necessary every year but with a low number of sessions.

6.4.11 Expected Treatment Outcome

Expected Treatment Outcome for Hair Removal

Fitzpatrick skin type	Hair Color	Fluence [J/cm ²]	Anatomic area of the body ²	Spot Size & Pulse Duration ¹	Average Number of Treatments ¹	Effectiveness (% of hair removal/regrowth)
I-III	Light brown to brown	33 - 44	Lower or upper legs	8 - 14 mm 2 x 60 msec	4 - 7	% of hair regrowth 67% with 8mm 54% with 12 mm 55% with 14mm
	Fair, brown to black	44	Lower or upper legs	8 mm	1 - 3	% of hair regrowth 11% ± 9% black hair 26% ± 22% brown hair 31% ± 28% fair hair

Fitzpatrick skin type	Hair Color	Fluence [J/cm ²]	Anatomic area of the body ¹	Spot Size & Pulse Duration ¹	Average Number of Treatments ¹	Effectiveness (% of hair removal/regrowth)
II-IV	Light brown to black	35	Upper lip	12 mm 2 x 45 msec	4 - 7	49% hair reduction
	Light to dark brown	25 - 35	Face, axilla, bikini	12 mm 50 msec	4 - 7	70 % hair removal
IV-VI	Dark brown to black hair	10 - 36	Plethora of body locations	12 mm	4 - 7	% of hair reduction 26% after 1 treatment 47% after 2 treatments 64% after 3 treatments

¹Individual response depends on hair type as well as biological factors that may affect hair growth patterns. Some patients may respond quicker or slower than the average number of treatments.

¹One cannot expect long term epilation in a single or even double treatment session. The duration of the resting period for hair follicles also depends on the body area.

Reference:

W. Baumlér, PhD; K. Scherer, MD; C. Abels, MD; S. Neff, MD; M. Landthaler, MD and R.-M. Szeimies, MD. The Effect of Different Spot Sizes on the Efficacy of Hair Removal Using a Long-Pulsed Diode Laser. *Dermatolog Surg* 2002; 28: 118-121

E. J. Fiskerstrand, MD, PhD; L.O. Svaasand, PhD; J.S. Nelson, MD, PhD. Hair Removal with Long Pulsed Diode Lasers: A Comparison Between Two Systems With Different Pulse Structures. *Lasers in Surgery and Medicine* 32: 399-404 (2003)

N.S. Sadick, MD, FACP, FAACS; V.G. Prieto, MD, PhD. The Use of a New Diode Laser for Hair Removal. *Dermatol Surg* 2003; 29: 30-34

Albert J. Nemeth MD, FAAD, FASLMS, FASDS, FISHRS. Impressive Laser Hair Removal with an 810 nm Powerpulsed Diode Laser: Safety and Efficacy. Presentation EADV 2004

6.5 Vascular treatment

A smaller spot adjusted to the vessels (optional handpiece) is used for the laser treatment of telangiectasias.

The laser treatment of superficial dilated vessels, particularly in the face but also in other locations, has been performed very successfully for more than 20 years – especially by using laser devices that emit light in the yellow and green ranges (argon, dye, KTP and copper bromide lasers).

These lasers are very well suited for the coagulation of very small vessels because their wavelength is absorbed very strongly by blood but very slightly by the water of the surrounding skin (see absorption curve). The selective photothermolysis method ensures an effective damage of the blood vessel and the simultaneous protection of the skin surrounding it. The penetration depth of this laser radiation is sufficient for very small vessels.

Due to the deeper penetration of laser light, lasers having wavelengths in the infrared range – like the diode laser The Epi Lab – are better suited for vessels with a larger diameter, for example for typical leg veins that are mostly larger and located deeper under the skin surface, and for hemangioma. The lower water absorption compared to the Nd:YAG laser (1064nm) is advantageous because normally fewer side effects develop.

Recent reports describe the successful treatment of port wine stains by using lasers in the near infrared range, which did not react well on the treatment with lasers in the green and yellow ranges before.

The relatively large spot diameter (that is possible due to the high power of The Epi Lab device) plays also an important role here: the penetration effectiveness of the laser radiation is the higher the larger the spot.

Leg vein telangiectasia are still a frequently asked cosmetic problem. The method of choice is sclerosing. The effectiveness of this therapy has been shown for decades. However, many patients prefer a laser treatment.

Generally, the success rate of the laser treatment of leg veins is lower than the one for treating, for example, telangiectasias in the face because they are under hydrostatic pressure and vary considerably in their diameter and oxygen saturation. The success rate that can be achieved is about 70%.

The absorption of laser light in melanin in dark-skinned patients is less for the wavelength of The Epi Lab than for lasers operating in the green-yellow range (see absorption curve). To minimize undesired side effects it must nevertheless be considered during the treatment by a sufficient cooling of the skin surface.

Apart from the wavelength, the energy density (FLUENCE) and the exposure time of the laser radiation (pulse duration) are the crucial parameters for a successful treatment. The fluence must be selected so that on the one hand the vessel will be closed and on the other hand undesired side effects are not caused (therapeutic window).

The pulse duration should be similar or shorter than the thermal relaxation time of the vessel (this is the time the target needs to emit 50% of the heat produced by the absorption) to follow the treatment concept of selective photo-thermolysis mentioned above. If the pulse duration is too long, the skin surrounding the vessel will be heated up more, the effectiveness will be reduced and the risk of side effects will increase.

6.5.1 Selection of patients

Patients with skin type I to III who are not tanned in the area to be treated can be treated with the following restrictions.

The Epi Lab is best suited for treating dilated vessels having a diameter of about 0.5. to 2 mm. The treatment of hemangiomas is also possible.

Leg veins can be differentiated into red and blue ones. The red vessels are smaller, more superficial and exhibit higher oxygen saturation (dilations of the arterial loop of the capillaries).

The blue ones are a dilatation of the venous capillaries with a lower flow speed. The larger vessels show a better reaction to the treatment than the smaller ones by means of this laser.

Before treating leg veins a varicosis, particularly of the lateral branch veins, is to be excluded, e.g. by Doppler sonography, or it is to be surgically treated or sclerosed.

Reticular varices should only be treated by laser if they occur solitarily or without connection to the depth

Reticular varices should only be treated by laser if they are identified as solitary varices or if they do not go deeply.

6.5.2 Contraindications

For contraindications including those for leg veins please see chapter *Intended Use*.

CAUTION

The higher fluence required for the treatment of dilated vessels and the stronger absorption due to the size of the vessel diameter cause more intensive side effects (e.g. compared to hair removal), subsequent pigmentary changes in particular.

Therefore, patients with skin type IV, V or VI must not be treated at all and tanned areas (postpone the treatment).

Useful information

Heavy smoking and/or overweight may adversely affect the efficacy of the treatment.

6.5.3 Possible side effects

The higher fluence required for the treatment of dilated vessels and the stronger absorption due to the size of the vessel diameter can cause more intensive side effects (e.g. compared to hair removal). Therefore, sufficient cooling of the skin as described in the following is absolutely necessary and more important than for hair removal treatments.

All side effects are less strong for patients with light skin and thinner vessels and become stronger the thicker the vessel and the darker the skin and, of course, if the fluence used is too high for the individual case (test treatment!).

Therefore, the start values for the test treatment of skin type III must be lower (darker skin types must not be treated).

During the treatment the patients have a pricking feeling connected with heat. It can even be a feeling of burning in case of thicker vessels.

After the treatment an erythema can be seen in the course of the vessels. If the vessels are thicker the formation of an oedema along the vessels is also possible. These effects are normal and last for a period of between few hours and 2 days. Blisters, incrustation and temporary pigmentary changes develop sometimes, but more often compared to epilation (light protection!).

Normally, blisters and incrustation heal without any negative consequences during about one week. Pigmentary changes can last for a period of up to some months.

Persistent pigmentary changes and scars are extremely seldom if the device has been operated as intended and if the treatment has been individually adjusted to the patient.

6.5.4 Preparing the laser treatment

As a compulsory pre-condition, the area to be treated must be cleaned. Remove also cosmetic products. The hair, if existing, must be carefully removed.

Photos of the initial situation should be taken to document the success of the treatment.

The treatment can cause pain, particularly for thicker vessels. However, anesthesia is normally not necessary if the area is sufficiently cooled before and after the treatment. The use of topical anesthetics, such as EMLA, does not reduce the sensation of pain.

A thin layer of a transparent, colorless ultrasound gel is to be applied onto the treatment area.

Cold cool packs (from the fridge, not frozen!) are to be laid out.

6.5.5 Test treatment

If the patient decides in favor of a laser treatment after getting the information from the physician, a test treatment should be performed afterwards in any case.

Due the high number of different vessel diameters the practically optimum fluence can only be found empirically by treating test areas and must be defined anew for each area that has another vessel diameter (color/structure change).

CAUTION

The higher fluence required for the treatment of dilated vessels and the stronger absorption due to the size of the vessel diameter can cause more intensive side effects (e.g. compared to hair removal).

Therefore, a sufficient cooling before and after the treatment is very important to minimize the side effects!

After preparing the treatment area, the test area to be treated is pre-cooled during some seconds by using a cool pack (ca. 4-8°C). Afterwards the cooling probe of the handpiece is moved backwards and forwards in longitudinal direction over the vessel to be treated without an active laser.

The test fluence for skin type III must be set lower than for light skin, the pulse duration as short as possible for the set fluence for smaller vessels and longer for larger vessels and skin type III. The treatment speed principally used at the beginning is set to 0.5 Hz (SPEED). This setting allows to assess the situation after only one pulse because it is possible to deactivate the foot switch after one pulse.

Moreover, the lower the treatment speed the longer the pre-cooling time for the next area to be treated.

WARNING

Laser radiation can cause serious damage to the eye.

Everyone present in the laser room during a treatment session must wear laser safety goggles according to the instructions given in chapter *Technical Data*!

The set fluence must be controlled before and during the treatment!

The Epi Lab is switched to the *READY* mode with the desired fluence.

Tip

Contrary to epilation treatments, pressure is **NOT** applied by the handpiece. .

It is recommended to keep a distance of about 1mm between laser probe and skin. At the same time the cold probe of the skin cooling has to be in contact to the skin (please tilt the handpiece a little bit).

2-3 pulses are applied to the area to be tested; the laser tip is moved towards the skin cooling. Afterwards, the area treated is assessed.

If the vessel has faded (more typical for smaller vessels) or becomes darker (more typical for larger vessels) the fluence will be kept and spot by spot is placed along the vessel.

Immediately after the laser treatment, the cooling probe is used to post-cool the area for some seconds.

If a reaction cannot be observed, wait some minutes to see if an inflammation reaction develops. In this time it is also possible to test another vessel size in the same procedure.

The larger the vessel the lower the fluence required normally.

If even after waiting a reaction cannot be seen, the NEW (!) test area is again pre-cooled by using the cooling probe and then it is lasered in an overlapping manner or the same area twice before the fluence is increased, if required. The risk of side effects is less in case of overlapping than for a too high fluence!

It is important that the test areas are pre- and postcooled again and again by the cooling probe.

After the treatment the laser is to be switched to the *STANDBY* mode.

By performing a test treatment it is also possible to assess the risk of pigment changes, curing and the risk of subsequent scarring.

Moreover, the patient can get an impression about the way in which the lesion will be treated.

6.5.6 Treatment procedure

The handpiece with integrated skin cooling is placed onto the skin so that the area to be treated first (already precooled by using the cool packs) is again cooled by the metal probe. The handling is very easy, because the skin cooling system is adapted to the laser handpiece.

The first area to be treated is pre-cooled for some seconds by using, for example, cool packs. Then, the first vessel is pre-cooled by the cooling probe as described above.

The first area is lasered with the parameters defined during the test (possibly different for different vessels) and then the probe is used for post-cooling the treated vessel.

If one area has been completely treated, a cool pack is applied for post-cooling. Then, the next area is treated in the same way.

It must be ensured that the probe is always kept in contact to the skin (without applying pressure) and that the handpiece is always moved towards the cooling probe.

After the treatment the laser is to be switched to the *STANDBY* mode.

The total treated area is continuously cooled for at least 10 minutes.

Due to the consequent cooling this procedure leads to a reduced sensation of pain and to a reduction of the side effects.

6.5.7 Post-treatment

After cooling the skin sufficiently, apply an anti-inflammatory ointment that should not be removed before next day.

In case of over-treatment a topic antibiotic is indicated and it is absolutely necessary to control the curing process.

CAUTION

If the treated area is exposed to sun and solar radiation after the treatment, stronger side effects can develop, pigmentary changes in particular.

Therefore, such a radiation should be avoided for a period of 6 weeks or a sun protection cream with a protection factor not below 20 should be used at least.

The laser area should not be traumatized either for 2 to 3 days.

Moreover, the patients must do without extreme physical or vessel stress, such as jogging, sports in the fitness studio, Finish baths, hot baths, etc. for a short time.

6.5.8 Follow-up treatments

If the vessels have not directly faded press your finger along the vessel and see whether blood is flowing back or not to assess the results.

Depending on the result of the treatment, the individual skin reaction and the interests of the patient, a further treatment with higher energy density can be scheduled, if required.

The interval between the sessions should not be shorter than 12 weeks. It is reported that even after some months vessels disappear. Up to this time, some vessels can be darker than before (blood coagulation). The patients must be informed about these facts.

More than three treatments are not useful for leg veins. As mentioned a success rate of 70% can be reached.

6.5.9 Expected Treatment Outcome

Study	Methods	Results
V. Tsepkenko, S. Khaustov, A. Burov "Our experience in vessels' laser coagulation" Clinic of plastic and aesthetic surgery "VIRTUS" Odessa, Ukraine (translation from Russian)	Objective: efficacy of selective laser coagulation in patients with blue veins (up to 3mm) Parameter: MeDioStar (Asclepion, Germany), wavelength 810 nm, spot size 6 mm, fluence 80-160 J/cm ² Patients: 106 patients with blue vessels up to 3mm diameter Patients of the period of last 6 months	complete healing – 75% (79 patients), significant improvement – 20% (21 patients), moderate improvement – 4% (4 patients), inconsiderable improvement – 1% (1 patient).
Wollina U, Konrad H, Schmidt WD, Haroske G, Astafeva LG, Fassler D. Response of spider leg veins to pulsed diode laser (810 nm): a clinical, histological and remission spectroscopy study. J Cosmet Laser Ther 2003;5:154-62	35 patients with spider leg veins, 2 treatments 2 weeks apart, 810 nm diode laser 12 mm spot size, 60 ms pulse length, 80-100 J/cm ² fluence Skin biopsies Remittance spectroscopy	Complete clearance after 1 st treatment with 15 of 35 patients, 20 of 35 remarkable improvement; stable situation after 6 month, 2 patients had mild scarring, otherwise no side effects.

6.6 Treatment of acne

Acne vulgaris is one of the most frequent diseases treated in dermatological practices. Mostly, the treatment takes a long time and is partly connected with considerable side-effects. Most of the young patients concerned suffer much from it.

In most of the cases it is a disease of the pilosebaceous follicles that are mainly located in the face, on the back and on the chest. First, they generate non-inflammatory blackheads but in the later stages a number of inflammatory efflorescences such as papules, pustules and nodules can also be developed. Acne vulgaris or 'usual acne' is caused by an increased androgen influence in the puberty and fades away until the beginning of the third life decade at the latest.

The infrared radiation of the EpiLab device penetrates – like necessary for hair removal – sufficiently deep into the skin to reach the pilosebaceous follicles.

The treatment as well as contraindications, side effects and treatment procedure are identical to epilation (see the corresponding chapters).

But unlike the epilation treatment, the acne treatment consists of 3-4 passes in direct succession and then it is the best to continuously cool the skin with a cooling mask. Due to the inflammation the treatment is painful.

3 to 5 treatments at one-week-intervals are required.

Acne inversa¹ is a disease of the skin appendage. It starts from the sebaceous glands and secondarily affects the apocrine perspiratory glands. Men and women are equally affected. This disease often starts during puberty but can also occur at any time later.

Predilection areas are especially axillas and the groin-pubic region, also spreading to the gluteal. The back and belly are less frequently affected and seldom the zones typical for acne like the face.

At the areas affected, first nodular, solitary foci of inflammation develop and in the course of the disease, they show recidivism and confluence to formations similar to carbuncles. Often, flat scar plates with a thickness of centimeters are caused. They have cutaneous/subcutaneous fistulas that can only be removed by an operation. The stress for the patient is not only caused by the inflammatory pains, the motoric restrictions by scar formation, but also by the repeated fetid infections.

The reasons for this disease may be a genetic disposition, its development in combination with other diseases (M. Dowling-Degos), keratinization defects of the follicular ducts, seborrhea, a special bacterial colonization by propioni bacteria and staph. Aureus, smoking and adiposity are also responsible for the development of this disease.

According to the medical history, in some patients the development of the disease can be possibly explained by an intensive hormone fluctuation (temporary secondary amenorrhea due to a considerable loss in weight), but the disease does not stop at the end of this possible cause.

The therapeutic measures taken so far are the repeated antibiotic and antiseptic treatment, abscess splitting as well as scar excision and defect coverage by means of split-skin graft. The disadvantages of this procedure are the repeated development of recurrences, the repeated small- or large-area surgical procedures, pain, the required bandage changes, the work time losses due to sick notes as well as the psychological stress with possible subsequently limited social contacts.

An alternative therapy option is the laser treatment by using a diode laser such as the EpiLab. The laser primary designed for epilation purposes penetrates deep under the skin surface and its effect is due to the absorption of light most of all by pigments or similar light-sensitive structures. The hair follicles are affected first. They are brought to involution due to the produced thermal damage after

¹ by courtesy of Dr. med. C. Passos Pereira, Dermatologe, Hurth-Hermuhlheim
published in DERMAforum, September 2003

several treatment sessions. Porpyrhin-forming bacteria and fungi, that are considered local causes of acne inversa are also destroyed in this procedure. After a normal epilation treatment, this effect becomes obvious by a reduced smell development of the axillas.

The treatment of the affected areas is performed at 810 nm with 30-40 J/cm². First, it is a single treatment of the affected area. In the most favorable case, the treatment should be carried out in the primary pustule or abscess stage without scar formation.

20 patients were treated within the last three years. Their diseases showed very different development stages starting from the recessive abscess formation in the inguinal zone up to the extended finding with partly considerable complications (repeated operations, condition after sepsis with vertebral abscess formation).

All patients showed a positive reaction to the therapy. At best, new abscesses healed within few days without requiring further antibiotic therapeutic measures. In case of inflammation reactions that have already existed for a longer period of time (longer than 6 months) the inflammation could heal within 6 weeks. In case of extended scar abscess cavities, an antiseptic effect could be achieved but a healing of the cavity not; recessive symptoms are to be expected here.

Patients who have not had extensive scars yet or areas not scared yet showed a good recessive prophylaxis. In the longest period about three years, a complete freedom of recessive symptoms could be observed after an abscess separation performed under general anesthetic once a month before. As an additional finding the smell development was reduced.

Compared with earlier investigations it can be said that a further therapeutic option for treating acne inversa exists. It can be just used at the beginning of the development of the disease and go without an operation in favorable cases.

For advanced disease cases, a significant improvement of the finding and the general condition of the patients can be achieved. Here, primary radical operations with defect coverage are not required and instead of this only minimum operations are necessary so that the disease can be better tolerated by the patients.

Expected Treatment Outcome

Study	Methods	Results
Shah GK. Efficacy of diode laser for treating acne keloidalis nuchae. Indian J Dermatol Venereol Leprol 2005; 71: 31-4	2 patients with acne keloidalis nuchae 810 nm diode laser 100 ms pulse length 23-26 J/cm ² fluence	Patients had about 90% clearance one month after 4 th treatment. Papules decreased in number and reduced in size. No new lesions appeared after 6 months. No side effects.
Passos Pereira C. Minimaleingriff ersetzt radikale Operation. Dermaforum 2003; Nr. 9, S. 22	20 patients over 3 years, Treatment analog epilation Minimum side effects	Regression of fresh inflammations within few days, of scary inflammation in 6 weeks

6.7 Skin rejuvenation

Generally, the term skin rejuvenation is understood as an improvement of the skin structure (lifting, reduction of the pore size), clearance of pigmentations and reduction of small wrinkles so that the overall appearance is rejuvenated.

Also here, the treatment, contraindications, side effects and treatment procedure are similar to epilation procedures. (See the corresponding chapters.)

Three treatments at 4 weeks intervals are to be performed.

But unlike the epilation treatment, the treatment consists of 3 passes in direct succession and then it is the best to continuously cool the skin with a cooling mask.

3 treatments at one-month-intervals are required.

If, however, age spots shall be brightened (diagnosis must be absolutely reliable to exclude cancerous and precancerous lesions!) a small spot size at a distance of about 1 mm from the skin is to be used for the treatment. After the treatments the lentigines become darker first! Make sure to inform your patients accordingly.

Lecture given by Sue McCoy, M.D. about Joint Laser Conference in Edinburgh, 2004

Translation from English:

Treatment of lentigines with a diode laser at 808nm

S. McCoy; Laser, Skin & Vein Clinic, North Adelaide, South Australia

Background

Lentigines are common benign epidermal skin lesions which occur predominantly on sun-exposed areas of individuals with fair (Fitzpatrick Types I-III) skin, increasing incidence with advancing age. Histologically, solar lentigines exhibit elongation of the epidermal rete ridges and an increase in the number of melanocytes and the amount of epidermal pigment.

Many different lasers have been used with a high degree of efficacy and low incidence of adverse effects including Q-switched Nd:YAG, ruby and alexandrite, frequency-doubled long-pulsed Nd:YAG, Erbium YAG, ultra-short pulsed CO₂ and others. This study evaluated treatment of lentigines with a long-pulsed 808 nm diode laser which is in common clinical use for the purpose of hair reduction.

Materials and method

Twelve healthy consenting female adults, age range 44 to 74 years, exhibiting numerous (at least 10) lentigines on the dorsum of the hands and/or the arms participated in this trial. A long pulsed 808 nm diode laser (The Epi Lab, Carl Zeiss Meditec, Jena, Germany) was used only in the 'Pro 2' mode – single pulse without interval. Three patients were initially treated with a 12 mm spot. The remaining 9 patients were treated with a 6 mm spot which allowed higher fluence with a shorter pulse width. Treatment was carried out with a silicon spacer to hold the circular flat lens 1 mm from the skin surface (in contrast to the recommended method of treatment for hair reduction which is skin contact with clear gel coupling).

All patients were treated with energies commencing at 25 J/cm² and increasing in 5 J/cm² increments until darkening of the lentigine within 5 minutes of treatment was observed. The remainder of the lesions were treated at this fluence. One lentigo per patient located on the posterior arm was treated at a pre-selected fluence and was biopsied immediately after the treatment. A total of 12 biopsies were studied with 2 at each fluence from 25 to 50 J/cm². At follow-up 4 weeks later the estimated percentage disappearance of lentigines was recorded as 0-25%, 51-75% and 76-100%. Any adverse effects were noted.

Results

All patients showed 50% or greater reduction in lentigines at one month except one patient treated with a 12mm spot and 30 J who showed only 25-50%. There was poor correlation between clearance and fluence but strong correlation between clearance and color of the original lesions with darker lesions showing greater clearance.

There was no scarring at one month but some lesions showed residual erythema at this short post-treatment interval. Histology consistently showed streaking of the nuclear chromatin of the pigmented basal keratinocytes, clefting of a thin upper rim of the papillary dermal collagen. The effects were more pronounced with the 6 mm than the 12 mm spot.

Summary

A diode laser emitting 25 to 50 J/cm² at 808 nm is capable of inducing limited and focal necrosis of solar lentigines in the forearm skin of Caucasian females with satisfactory clinical results. Treatment is rapid and well tolerated. Results are more dependent on the color of the lentigino than the fluence of the laser. Operator judgement and experience is therefore important in obtaining good results.

Expected Treatment Outcome

Study	Methods	Results
Kilmer S, Chotzen VA, McClaren M. Use of the LightSheer™ Diode Laser System for the Treatment of Benign Pigmented Lesions. ©2000 by Lumenis Inc	24 patients with 69 pigmented lesions, skin types I-IV, 800 nm diode laser, 9x9mm spot size 15-30 ms pulse length, (double pulse) 30-40 J/cm ² fluence in the beginning (later up to 60 J/cm ²)	26% of patients had >75% improvement after 1 st treatment, 67% had <50% improvement, One month after 2 nd treatment 57% of patients had >75% improvement, 32% had <50% improvement, Some patients had temporary side effects (pigmentation changes)
Dr. Dai, Dr. Lai "Clinical evaluation of skin rejuvenation on Chinese skin with MeDioStar" Taipei, Taiwan	Objective: Evaluation of the effect of MeDioStar treatments for skin rejuvenation on Chinese skin Parameter: MeDioStar miXT, spot size 12mm Wavelength 810 / 940nm, Fluence 30 – 40J/cm ² Pulse duration 300 – 400ms 3 passes, 1 treatment Evaluation: Objective measurement of skin of 6 patients concerning: Pore size, wrinkles, texture, brown spots, by a special camera (Canfield Visia System) with objective measure	patient 1: (2 weeks after 1 st treatment) skin texture improvement: 54 to 85 pore size improvement 45 to 60 patient 2: spot (pigment) improvement 24 to 39 patient 3: (immediately after 1st treatment) pore size improvement 32 to 51 wrinkle improvement 60 to 99 skin texture improvement 29 to 51 patient 4: (1 day after 1st treatment) pore size improvement 3 to 13 patient 5: (15 days after 1 st treatment) pore size improvement (64 to 79) patient 6: (6 days after 1st treatment) wrinkle improvement 66 to 99


7 Accessories

WARNING

The use of accessories not tested and released by the manufacturer for the use with this laser device can cause severe damage to the eye or electrical shocks.

Only use accessories and cables provided by Asclepion Laser Technologies for this device.


Important information

We strongly advise against the use of accessories made by other manufacturers. Even if an official testing authority should have certified that a specific accessory unit could be used safely, Asclepion Laser Technologies cannot assume any liability for these products.

On request, Asclepion Laser Technologies will readily send you an updated list of accessories.

8 Cleaning, disinfection and sterilization

Important information

The following procedures for cleaning, disinfection and sterilization must be observed absolutely to keep the device and accessories in best shape. Any deviation from the specified procedures is prohibited and of course any opening of the device or handpieces - otherwise any warranty expires.

8.1 General rules

Device

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock!

Detach line power cable from the mains socket in order to break power supply to the system before you proceed with work for cleaning and disinfection.

Only outer system surfaces may be cleaned by the user. For cleaning and disinfection solely damp wipes have to be used, e.g. Kodan® wipes. The cleaning or disinfection using wet clothes or the spraying e.g. with disinfectants is prohibited.

Handpieces

WARNING

The contact surface of the handpiece contacts the patient's skin during the treatment. There is the risk of microorganisms being transferred.

Handpiece parts must be cleaned and disinfected immediately following use and before application on a patient.

If it should be necessary to clean the device during a treatment the device must be switched to the STANDBY mode before starting cleaning works! Do not use alcohol or isopropyl alcohol near the point of treatment because this causes risk of fire if the laser is activated afterwards!

CAUTION

Dirt settling on an optical surface may burn in when exposed to laser radiation and cause side effects such as reddening, possibly blistering, incrustation, seldom scars!

Ensure that all accessible optical surfaces are not damaged and clean after the cleaning process and remain clean during the treatment.

The procedure instructions given in the following sections are to be followed.

8.2 General procedure for cleaning and disinfecting the handpieces

Important information

You are prohibited from cleaning a handpiece with wet clothes or under running water or from spraying e.g. with disinfectants. This may destroy the internal optics and electronics. Solely moistened clothes have to be used - for disinfection e.g. Kodan wipes. (The parts are damaged by strongly alkaline solutions!)

The handpieces cannot be sterilized.

Frequent cleaning and disinfection have low influence on the accessories. The color may become lighter or changes but this does not affect the use.

The end of the product service life is normally determined by wear, tear and damage caused by its use.

If the handpieces are not in use for a longer period they must be protected against contamination.

INSTRUCTIONS

Workplace	The workplace must be clean and disinfected.
Preparations for cleaning	The cleaning of the accessories should be performed as soon as possible after their use. For this purpose they are to be removed from the device. Opening or dismantling of the handpieces is prohibited.
Cleaning: manual	Equipment: microfiber wipe delivered by Asclepion or similar lens cleaning tissues; disinfection wipes (e.g. Kodan® wipes). When using the wipes follow the manufacturer's instructions 1. Procedure for the body of the handpiece: The accessible surfaces of the handpiece body are cleaned by using a soft, slightly damp cloth or by the disinfection wipe. The cloth must not be wet or even dripping. You can use neutral cleaning agent for stubborn dirt. 2. Procedure for optical surfaces: The accessible optical surfaces must be absolutely clean and lint-free. Clean them only, if this is not the case. For doing this, first the surface is breathed on and then cleaned with a microfiber (as indicated above). For persistent dirt disinfection wipes have to be used. The optical surface must not get scratches and afterwards it must be lint-free and free of streaks.
Cleaning: automatic	not possible
Disinfection	The accessible outer surfaces are to be disinfected by means of a disinfection <u>wipe</u> (e.g. Kodan cloth). The optical surface must not get scratches and afterwards it must be lint-free and free of streaks
Sterilization	not possible, not required

INSTRUCTIONS

Storage and transport

The parts must be stored in a clean condition and protected against contamination. Ensure that the parameters of the correspondingly admissible ambient conditions indicated in chapter *Technical Data* are given.

The manufacturer has validated that the above instructions are **SUITED** for preparing the accessory parts for their reuse.

WARNING

The contact surface of the handpiece contacts the patient's skin during the treatment. There is the risk of microorganisms being transferred.

The person who cares about cleaning has all the responsibility that the actually performed cleaning and disinfecting works by means of the equipment and materials and with the staff involved achieve the desired results.

To guarantee this, validations and routine monitoring of the procedures are normally required. The effect and adverse consequences of each deviation from the instructions given should also be evaluated.

8.3 Special procedure for cleaning and disinfecting the handpieces


Figure 34 Handpiece

During the treatment:

The optical contact surfaces of the laser tip are largely protected against contamination to a great extent by applying the recommended technique. Nevertheless they are to be kept in clean condition during the treatment, too (particularly when treating thicker dark hair)!

CAUTION


Dirt settling on an optical surface may burn in when exposed to laser radiation and cause side effects such as reddening, possibly blistering, incrustation, seldom scars!
Ensure that the laser tip of the handpiece remains clean!

For this reason, the optical contact surface must be checked regularly and cleaned as instructed above. You can, for example, rapidly wipe off the surface each time when a line has been lasered before starting to laser the next line. This method avoids that a possibly existing contamination burns in when the laser treatment is continued.

After the treatment:

The device must be switched off and the handpieces have to be cleaned and disinfected as described above immediately thereby the dirt can be removed easily..

9 Fault messages and fault removal

Important information

The following procedures for fault removal have to be observed absolutely. Any deviation from the procedures described herein is prohibited and of course any opening of the device or handpieces – otherwise any warranty expires.

9.1 General information

If the device does not respond to turning the key switch to ON, make sure that the following requirements are met:

1. The device is connected to the line power supply.
2. The power switch at the rear panel of the device is switched on (position 'I', see *Figure: Connections at the rear panel of the device*).
3. The Laser emergency STOP switch is unlocked. If locked, please unlock by turning the red knob anticlockwise until the lock is released and the switch jumps out (see *Figure 5*).

As mentioned before, you are assisted in fault finding by operator guidance via the display.


Figure 35: Fault messages screen (example)

For general faults, the display shows a fault message with instructions in plain text. Some faults can be removed by the user. For doing this, follow the instructions on the display.

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock and damage to the eye!

You are strongly advised against trying to do repair works at the device yourself that go beyond the instructions given on the display!

In the case of safety-relevant faults, the device will shut down within milliseconds. A display of a fault message is impossible. If this happens or any other fault appears that has not been described in this manual, disconnect the device from line power by unplugging the power plug from the power outlet.

This chapter provides information on possible causes of any malfunction of The Epi Lab system. Faults can be removed by the operating personnel except of those specially referred to as absolutely requiring repair by technicians employed with or authorized by Asclepion Laser Technologies.

Improper use or maintenance of the laser system may void the warranty granted. Before starting any repair or if you have questions that have not been answered by this manual, contact your service partner (see chapter 10 *Customer service*).

9.2 Fault messages and fault removal

The Epi Lab is provided with alarm systems. If they are released, the system will switch to the safe mode and stop the emission of laser radiation. Please try to remove the fault yourself by following the instructions given in the table below. If the error persists, contact Asclepion Laser Technologies or a service partner expressly authorized by Asclepion Laser Technologies.

Before you contact the Technical Service, please make a note of the fault message and the error code displayed.

Fault message	Explanation / Corrective action
Communication error (all 50 XXX errors)	The communication is interrupted. Start the device anew. Operation can be continued. Please contact the Technical Service if the error occurs again.
Temperature error (all 150 XXX errors)	The operating temperature is too low or too high. Check if the device is sufficiently ventilated. Start the device anew. Operation can be continued. Please contact the Technical Service if the error occurs again.
Foot switch (400 XXX errors)	Please release the foot switch. (do not activate the foot switch too early, wait until the display shows READY with yellow warning sign)
Door interlock (600 601 errors)	The door interlock contacts are open or the door interlock plug is not connected. If the door has been connected to the door interlock, door open. Connect the door interlock contacts correctly or activate the door interlock switch. If the door has been connected to the door interlock, close the door.
Handpiece error (all 110X errors)	Check the handpiece connection. Start the device anew. Please contact the Technical Service if the error occurs again.

Possible malfunctions

Malfunction	Explanation / corrective action
No power supply	Verify that the power cable has been connected to a power outlet. Is the power switch set to position 'I'? Verify that the Laser emergency Stop switch is released. Check the line fuses of the domestic power supply system.
Function keys of the display do not react (or only sporadically)	Faulty calibration of the touchscreen The calibration is to be performed by an authorized service engineer. Therefore, inform your vendor or directly the Technical Service of Asclepion Laser Technologies (see chapter 10 <i>Customer service</i>).
No laser emission with the foot switch pressed	Verify that the connecting cable of the foot switch is properly connected. The device is in STANDBY. Press the READY key to switch to the READY mode.

10 Customer service

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock and damage to the eye!

Never open the device and attempt to repair the device yourself if the laser does not work properly.

Only specialists appropriately qualified and authorized are allowed to open the laser device.

The laser does not contain any user-serviceable components. Only service engineers employed with or trained by Asclepion Laser Technologies are authorized to service the device. This also applies to the prescribed annual safety tests.

If you need assistance by the Technical Service, please contact us via the address or phone number given below. When contacting your service partner, always keep the model and the serial number of the device at hand. These numbers are specified on the type plate located at the rear panel of the laser device.

Asclepion Laser Technologies GmbH

Service

Brüsseler Str. 10

D - 07747 JENA

Germany

Tel.: +49 (0) 3641 / 7700 - 401


Fax.: +49 (0) 3641 / 7700 - 402

E-mail: service@asclepion.com

Please note that the manufacturer, installer and importer will not consider themselves responsible for the effects of servicing on safety, reliability and performance of the device unless the following requirements are met:

- installation, extensions, readjustments, modifications or repairs have been performed by authorized persons,
- the electric installation conditions of the laser room meet the requirements of IEC 64, as amended, and
- the device has been used in accordance with the instructions given in this manual.

11 Regular maintenance, safety checks and calibration

Important information

The following procedures for maintenance, safety checks and calibration have to be observed absolutely. Any deviation from the procedures described herein is prohibited and of course any opening of the device or handpieces – otherwise any warranty expires.

The laser does not contain any user-serviceable components.

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock and damage to the eye!

The operation, maintenance, test and calibration instructions given in this manually are to be followed.


Only service engineers employed with or trained by Asclepion Laser Technologies are authorized to service the device. This also applies to the prescribed annual safety tests, calibration included.

CAUTION

A too high output power of the laser can cause side effects such as reddening, possibly blistering, incrustation, seldom scars!

The output power of the therapeutic laser must be checked by a service technician of Asclepion Laser Technologies or a qualified service technician of the distributor within the technical safety check at least once a year, the device is to be calibrated too if required. The results of the check are to be documented in the Medical Device Logbook.

Important information

The annual maintenance is mandatory for the extended warranty.

To achieve a high reliability and life time it is strongly recommended to use the cooling liquid specially developed by Asclepion for this device (corrosion protection, antibacterial, frost protection, effective cooling)!

This coolant should **never** be mixed with distilled water – please observe the relevant label at the device (please see chapter *Installation / Labels at the device*).

An actual list of consumables is available on request from customer service.

11.1 Routine maintenance

The routine maintenance procedures described below can be performed by the personnel of the clinic or hospital.

General

Inspection of the outer surfaces of the device

Check if all cables are intact and connected firmly.

Frequency: weekly

To be checked by: _____ personnel of clinic or hospital

Check of the cooling system

Check if all ventilation grids of the cooling system are free.

Frequency: weekly


To be checked by: _____ personnel of clinic or hospital

Check of the handpiece

Check if the handpiece is in perfect state.

Frequency: before each use

To be checked by: operator of the device

Important information

Before the real counter reaches end of life time (and "Remain" reaches 0) there are several reminders on the display that the attached handpiece has to be send in for revision to the service department of Your distributor in a foreseeable time.

This way a long term planning is possible. During revision it is decided, if the handpiece can be used further for a specific number of pulses or if the handpiece has to be replaced.

11.2 Regular safety check

SCOPE OF ANNUAL SAFETY CHECK

The safety check should be carried out in compliance with the relevant national legally regulations, as amended, and include at least the tests mentioned below. The tests should be run in the order specified:

- visual inspection of laser device and accessories
- measurement of protective ground wire resistance
- earth leakage test under normal conditions
- measurement of the actual output power of the laser modules at the handpiece; calibration of laser power, if necessary
- functional test.

Tip

To inform yourself about the extent and performance of the safety checks, please contact your distributor or directly Asclepion Laser Technologies.

Additional information


Upon request, the manufacturer will provide circuit diagrams, component parts lists, descriptions, calibration instructions, or other information not yet contained in this User Manual to assist appropriately qualified technical staff to repair those parts of the device that have been designated by the manufacturer as repairable. In this context, "appropriately qualified technical staff" means staff that attended a manufacturer's service training course on this device and that was authorized to repair it.

12 Disposal

The device is to be disposed of according to regulation 2002/96/EC of the European Council for Waste Electric and Electronic Equipment [WEEE] or the specific national regulations on waste electric and electronic devices, as amended.

Please consult our Technical Service if you have any questions regarding disposal.


A Appendix: EMC manufacturer's declaration

Important information


The usage of others than the following accessories, transducers and cables may result in an increased interference emission and/or decreased immunity of The Epi Lab.

Accessory / designation		Length / dimensions
Foot switch with cable	To be ordered from Asclepion	$l < 3.0 \text{ m}$
Cold device supply line	Laser Technologies or the authorized distributor	$l < 3.0 \text{ m}$

Important information

This device is subject to special terms of network connections. The maximally allowed impedance of the power grid is $Z_{\max} = 0.1753 \text{ Ohm}$ at the power supply point. If necessary, the local energy supplier is to be contacted to ensure that the device will be plugged in a supply which has an impedance lower or equal $Z_{\max}$.

Important information

The Epi Lab must not be stacked or arranged directly next to or with other devices. If The Epi Lab is used close to other devices it must be monitored to control its proper operation.

A.1 Electromagnetic emissions

The Epi Lab is suitable for use in the specified electromagnetic environment. The customer and/or the user of The Epi Lab should ensure that it is used in an electromagnetic environment as described below.

Emissions measurement	Compliance	Electromagnetic environment guidelines
High-frequent interference emission to CISPR 11	Group 1	The Epi Lab uses HF energy only for its internal function. Therefore, its emissions of HF interferences are very low and are not likely to cause any interference in nearby electronic equipment.
High-frequent interference emission according to CISPR 11	Class B	The Epi Lab is suitable for use in all establishments, including domestic establishments and those directly connected to the public power supply network that supplies buildings used for domestic purposes.
Emitted interference of harmonics according to IEC 61000-3-2	Class A	For harmonics: The Epi Lab is intended for use in all establishments, including domestic establishments and those directly connected to the public power supply network that also supplies buildings used for domestic purposes
Emitted interference of voltage fluctuations / flickers according to IEC 61000-3-3	Fulfilled 1*)	For flickers: During operation an influence of the supply network cannot be absolutely excluded. If necessary, it is to be ensured – after consultation with the energy supplier- that the max. admissible grid impedance $Z_{max} = 0.1753 \text{ Ohm}$ is not exceeded.

1*) acc. to IEC 61000-3-11:2000 with max. allowed impedance of the power grid of $Z_{max} = 0.1753 \text{ Ohm}$.

A.2 Electromagnetic immunity

The Epi Lab is suitable for use in the specified electromagnetic environment.

The customer and/or the user of The Epi Lab should ensure that it is used in an electromagnetic environment as described below.

Immunity test	Test level acc. IEC 60601-1-2	Compliance level	Electromagnetic environment guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	+/- 6 kV contact charge +/- 8 kV air discharge	+/- 6kV +/- 8kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Fast transient electrical interferences/bursts acc. to IEC 61000-4-4	+/- 2 kV for power lines +/- 1 kV for input/output lines	+/- 2kV *	Mains power quality should be that of a typical commercial and/or hospital environment.
Surges acc. to IEC 61000-4-5	+/- 1 kV series-mode voltage +/- 2 kV common mode voltage	+/- 1kV +/- 2kV	Mains power quality should be that of a typical commercial and/or hospital environment.
Voltage dips, short interruptions and voltage fluctuations of power supply input lines acc. IEC 61000-4-11	< 5 % UT for ½ cycle (> 95 % interruption) 40 % UT for 5 cycles (60 % interruption) 70 % UT for 25 cycles (30 % interruption) < 5 % UT for 5 s (> 95 % interruption)	< 5 % 40% 70% < 5 %	Mains power quality should be that of a typical commercial and/or hospital environment. If the user of the device requires continued operation during power mains interruptions, too, it is recommended to power the device by an uninterruptible power supply or a battery.
Magnetic field for the supply voltages (50/60Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	The magnetic field strength should be at levels typical for commercial and/or hospital environment.

Note: U_T is the AC mains voltage prior to application of the test level

* not tested, because line length shorter than 3m

The Epi Lab is suitable for use in the specified electromagnetic environment. The customer and/or the user of The Epi Lab should ensure that it is used in an electromagnetic environment as described below.

Immunity test	Test level acc. IEC 60601-1-2	Compliance level	Electromagnetic environment guidelines
Conducted HF interferences acc. to IEC 61000-4-3	3 Veff 150 kHz to 80 MHz	3V	Portable and mobile radio equipment, including cables, should be used not closer to any part of the device than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter Recommended separation distance: $d = \left[\frac{3,5}{3V} \right] \sqrt{P}$
Radiated HF interferences acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{3,5}{3} \right] \sqrt{P}$ for 80 MHz to 800 MHz $d = \left[\frac{2}{3} \right] \sqrt{P}$ for 800 MHz to 2,5 GHz

where P is the nominal power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^a

Field strength of fixed radio transmitters, as determined by an electromagnetic site^a survey should be less than the compliance level in each frequency range^b

Interference may occur in the vicinity of equipment marked with the following symbol:

(((•)))

^a Field strengths of fixed transmitters, such as base stations for radio telephones (cellular/cordless) and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance levels given above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b Over the frequency range from 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

A.3 Recommended protection distances

Recommended protection distances between portable and mobile HF telecommunications equipment and The Epi Lab according to IEC 60 601-1-2.

Transmitter frequency	150 kHz to 80 MHz	150 kHz to 800 MHz	800 MHz to 2.5 GHz
Equation	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
Rated output of the transmitter (Watt)	Protection distance (meter)	Protection distance (meter)	Protection distance (meter)
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the protection distance can be estimated by using the equation in the corresponding column, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

Note

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D

Identification number: MMo/13/05/07-06

Member code/year/month/day-sequential number

Page 1 of 10

FB 7302 RAaC, 60601-1-4 MeDioStar
1550 Vs 06 MMo/130507-06.doc

Medical Product: MeDioStar

Item number: 1590 - 1599

Assessed construction level: Serial Device

Revision user's Manual:

☐ MDD 93/42/EEC☐ EN 60601-1 2006☐ EN 60601-1-2 2001+ A1:2006☐ EN 60601-1-2 2007☒ EN 60601-1-4 1996+ A1:1999☒ EN 62304 2006☐ EN 60601-1-6 2004☐ EN 60601-1-6:2007☐ EN 62366 2008☐ before R&D☐ during R&D☒ after sales☐ review reg. changed standards

Related documents:

Drawn up (history next page):

Vs.	Description of change	Date	Name	Sign
06	Review	May 07, 2013	M. Moller	

Checked and released:

	Division	Date	Name	Sign
Participant 1	QS R&D	May 07, 2013	P. Moller	
Participant 2				
Released by	R&D	May 07, 2013	L. Streit	

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
Form 7302	Risk analysis and control acc. to EN ISO 14971:2009	 Page 2 of 10 FB 7302 RAaC 60607 1.4 MoD oStar 1590 Vs 05 MMo130507-06 doc
Revision: D	Identification number: MMo/13/05/07-06 Member code/year/month/day-sequential number	

Drawn up:

Vs.	Description of change	Date	Name	Sign
06	Review	May. 07. 2013	M. Möller	sgd. M. Möller
05	Review (especially for cooling hand piece)	Nov 15, 2011	M. Möller	sgd. M. Möller
04	Review (especially for non-existing data, wrong data, temporal data sequences)	April 13, 2011	R. Herbst	sgd. R. Herbst
03	Review	Feb. 09. 2011	R. Herbst	sgd. R. Herbst
02	Update risks, Add on risk control measures	Oct, 04 2010	M. Kleiner	sgd. M. Kleiner
01	First edition	June, 16 2010	M. Kleiner	sgd. M. Kleiner

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision D

Identification number: MMo/13/05/07-06

Member code/year/month/day-sequential number

Page 3 of 11

FE 1302 RANC 60601-1-4 MMD05127
1530 Vs 06 MMD 130507 29 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
1.5 No ionising radiation – laser radiation	1. Laser	Automatic release	Eye damage / injury, Injuries by fright	Critical	Occasional	ALARP	- Electronic shutter (footswitch controlled, software independent) - Electrical and EMC-safe construction of device	Critical	inconceivable	passable	Software verification of control system Test report DIN EN 60601-1-2
	2. Laser	Footswitch released unintentionally (laser radiation)	Eye damage / injury, Injuries by fright, Injury of third persons	Critical	Occasional	ALARP	- Standby-Ready circuitry - Optical and acoustic laser emissions indication - Deactivation of laser in periods of inactivity - Footswitch with cover	Critical	inconceivable	passable	Software verification of control system Verification of the device
	3. Laser	No interruption of laser radiation in case of footswitch release (in terms of footswitch malfunction and/or short circuit in signals)	Eye damage / injury, Injuries by fright, Injury of third persons	Critical	Occasional	ALARP	- Electronic shutter - Footswitch with two independent contacts - Foot switch inquiry before - Laser-Emergency STOP button	Critical	Hardly conceivable	ALARP	Software verification of control system Verification of laser

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision D


Identification number: MMo/13/05/07-06

Member code/year/month/day-sequential number

Page 4 of 11

FS 7302 RAaC 60601-1-4 MeDiCSia-
1590 Vs 06 MMo/13/05/07-06 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
	4 Laser	Unintended release of laser radiation at Ready activation because of active foot-switch	Eye damage / injury, Injuries by fright, Injury of third persons	Critical	Occasional	At ARP	- Foot switch inquiry before laser ready activates control signals for laser firing - Electronic shutter - Laser-Emergency STOP button	Critical	inconceivable	passable	Software verification of control system Verification of laser
	5 Laser	People in the room without safety goggles during laser or pulse light treatment	Eye damage / injury, Injury of third persons	Critical	Occasional	ALARP	- optical and acoustic display of emission 2 seconds lid closing time at Ready activation	Critical	Hardly conceivable	ALARP	Software verification of control system
	6 Laser	Access of third persons during laser treatment	Eye damage / injury, Injury of third persons	Critical	Occasional	ALARP	- Door Interlock - external warning lamp	Critical	inconceivable	passable	Software verification of control system
	7 Laser	wrong laser parameter through faulty software and/or hardware (abnormal system end)	- Over treatment - Vesicles and / or scars	Critical	Occasional	ALARP	- Watch Dog System - Gate locking - Laser-Emergency STOP button - Approved CPU with well-known operating system and no second supplier for CPU and operating system (set of drawings)	Critical	inconceivable	passable	Software verification of control system Verification of laser Verification of the device

Form 7302	Risk analysis and control acc. to EN ISO 14971:2009	
Revision: D	Identification number: MMo/13/05/07-06 Member code/year/month/day-sequential number	Page 5 of 11 FB 7302 RAaC 60601-1-4 MedDrSta 1590 Vs 06 MMo130507-06.doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
	8. Laser	unintentional switch-on of the laser system	• Eye damage • Over treatment • Vesicles and / or scars	Critical	Occasional	ALARP	• Laser STOP button • Main switch for current supply • Key switch • Standby at Power on	Critical	inconceivable	passable	Verification of laser Software verification of control system
	9. Laser	wrong laser parameter through wrong read or non existing data of non-volatile storage	• Over treatment • Vesicles and / or scars	Critical	Occasional	ALARP	• CRC by EEPROM data • Dual data management	Critical	inconceivable	passable	Software verification of control system
	10. Laser	wrong laser output parameter after change of the diode driver	• Eye damage • Over treatment • Vesicles and / or scars	Critical	Occasional	ALARP	• monitoring of drive value and pulse width • Efficiency of laser headpiece non-volatile stored • Efficiency of diode driver non-volatile stored • Efficiency adjustment after change	Critical	inconceivable	passable	Software verification of control system Instruction at the service manual, service training

Form 7302	Risk analysis and control acc. to EN ISO 14971:2009	
Revision: D	Identification number: MMo/13/05/07-06 Member code/year/month/day-sequential number	Page 6 of 11 FB 7302 RAaC 6/06C 1-4 MeDisStar 1590 Vs 06 MMo/13/05/07-06 600

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
	11. Laser	wrong laser parameter after change of the laser hand piece	• Eye damage • Over treatment • Vesicles and / or scars	Critical	Occasional	ALARP	• monitoring of drive value and pulse width • Efficiency of laser headpiece non-volatile stored • Efficiency of diode driver non-volatile stored • Efficiency adjustment after change	Critical	inconceivable	passable	Software verification of control system Instruction at the service manual, service training
	12. Laser	Malfunction of diode driver wrong output parameter (incorrect current or pulse width)	• Eye damage • Over treatment • Vesicles and / or scars	Critical	Occasional	ALARP	• CRC by EEPROM data • Dual data management • Gate locking	Critical	inconceivable	passable	Software verification of control system
	13. Laser	Undefined operating status through I ² C BUS malfunction	Eye damage	Critical	Occasional	ALARP	• I²C BUS protocol with malfunction message	Critical	inconceivable	passable	Software verification of control system

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision D

Identification number: MMo/13/05/07-06

Member code/year/month/day-sequential number

Page 7 of 11

FB 7302 RAaC 60607 1 4 MeDiStor
1590 Vs 06 MMs/13/05/07-06.doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
3.4 Insufficient coolant supply or limitation of cooling	1. Cooling hand piece	- Malfunction of two point regulation of the peltier - Overheating - wrong read or non existing data of temperature regulation	- Over treatment - Vesicles and / or scars	Minor	Occasional	ALARP	- permanent check of the water temperature of the hand piece cooling unit - permanent check of the temperatures of cold side at the cooling hand piece - permanent check of the water flow of the hand piece cooling unit - Check sum over files with setup data	Minor	inconceivable	passable	Software verification of control system
	2. Cooling hand piece	- Malfunction of two point regulation of the peltier - Under cooling - wrong read or non existing data of temperature regulation	edema at touching	Minor	Occasional	ALARP	- permanent check of the water temperature of the hand piece cooling unit - permanent check of the temperatures of cold side at the cooling hand piece - permanent check of the water flow of the hand piece cooling unit - Check sum over files with setup data - Limited cooling power of the peltier	Minor	inconceivable	passable	Software verification of control system Verification temperatures

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D


Identification number: MMo/13/05/07-06

Member code/year/month/day-sequential number

Page 8 of 11

F0 7302 RAaC 60601-1-4 MMedStar
1560 V5 06 MMo/13/05/07-06 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
3.5 Storage/ operating outside required environmental conditions	1. Laser	Heat sink temperature too high	- overheating of the device - deviation of laser parameters - over treatment	Minor	Occasional	ALARP	- permanent check of main water temperature for the diode and cooling liquid temperature of the compressor - permanent adjust of fan speed - advices in instruction manual	Minor	inconceivable	passable	Software verification of control system check list EN 60601-1
	2. device environment too wet	short circuits in the device caused by leakage connections	- risk of fire - risk of electric shock	Critical	Occasional	ALARP	- permanent check of the water temperature of the hand piece cooling unit - permanent check of the temperatures of cold side at the cooling hand piece - permanent check of the water flow of the hand piece cooling unit	Minor	inconceivable	passable	Software verification of control system check list EN 60601-1
3.6 Incompatibility with combinations of other devices	n/a	n/a	n/a	n/a	n/a	n/a	There is no action needed, because it's no combination provided	n/a	n/a	n/a	n/a

Form 7302	Risk analysis and control acc. to EN ISO 14671:2009	
Revision: D	Identification number: MMo/13/05/07-06 Member code/year/month/day-sequential number	Page 9 of 11 FB 7302 RAaC 60601-1-4 MeOoStar 1590 vs 06 MMo130507-06.doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
5.11 Erroneous measurement and other measurement aspects	1 Laser	wrong laser parameter by wrong hand piece	- Over treatment - Vesicles and / or scars	Minor	Occasional	ALARP	- Hand piece data, saved in EEPROM with hand piece-type display and/or diameter, shown at the display - Protection against hazardous initial value through displayed energy density is shown at the display	Minor	inconceivable	passable	Software verification of GUI
	2 Laser	wrong laser parameter through wrong current provoked by malfunctions in diode driver or control value	- Over treatment - Vesicles and / or scars	Minor	Occasional	ALARP	- monitoring of drive value and pulse width - gate locking - redundant current measurement	Minor	Hardly conceivable	passable	Software verification of control system
	3 Laser	wrong laser parameters because of wrong display	- Over treatment - Vesicles and / or scars	Minor	Occasional	ALARP	Display communication with confirmation	Minor	inconceivable	passable	Software verification of control system
	4 Laser	wrong temporal sequence of master-slave communication	- Over treatment - Vesicles and / or scars	Minor	Occasional	ALARP	Display communication with confirmation	Minor	inconceivable	passable	Software verification of control system

Form 7302	Risk analysis and control acc. to EN ISO 14971:2009	
Revision: D	Identification number: MMo/13/05/07-06 Member code/year/month/day-sequential number	Page 10 of 11 FR 7302 RAx C 60501-1-4 MeDiGStar 1590 vs 06 MMo/130507-06 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
	5. Laser	wrong temporal sequence of laser pulse and evaluation	- Over treatment - Vesicles and / or scars	Minor	Occasional	ALARP	Synchronized laser pulse with evaluation and monitoring of overrun	Minor	inconceivable	passable	Software verification of control system
	6. Laser	wrong laser parameters because of wrong pulse length	- Over treatment - Vesicles and / or scars - Third-degree burns on to skin	Minor	Occasional	ALARP	- monitoring of the clock frequency - gate locking	Minor	inconceivable	passable	Software verification of control system
7.9 Failure of the device/function because of function/assembly faults	1. Laser	failures of current sensors, failures of skin cooling	- Overtreatment - Vesicles and / or scars	Minor	Occasional	ALARP	- Permanent monitoring of hardware status and software status, watchdog system - redundant current measurement - Automatic device standby / turnoff in case of failures - Laser-Emergency button	Minor	Hardly conceivable	passable	Software verification of control system
8.1 Hazards by investigation in the FDA database	1. Laser	energy output too high	- Overtreatment - Vesicles and / or scars	Minor	Occasional	ALARP	- redundant current measurement system - gate locking - annual control of laser output calibration	Minor	Hardly conceivable	passable	Software verification of control system User's manual

Form **7302****Risk analysis and control acc. to EN ISO 14971:2009**

Revision: D

Identification number: **ZI/11/02/03-01**

Member code/year/month/day-sequential number

Page 1 of 7

FB 7302 Risk Analysis and control
62601-1-6 MDS SE Vs 93 docD

Medical Product: Mediostar (NeXT)

Item number: 1590 ff

Assessed construction level: Pre-Serial Device

Revision user's Manual: 1590 99 001 Vs.05
Jan 11, 2011☐ MDD 93/42/EEC☐ EN 60601-1:2006☐ EN 60601-1-2 2001+ A1 2006☐ EN 60601-1-2 2007☐ EN 60601-1-4:1996+ A1 1999☐ EN 62304:2006☐ EN 60601-1-6 2004☒ EN 60601-1-6 2007☒ EN 62366:2008☐ before R&D☒ during R&D☐ after sales☐ review reg changed standards**Related documents:****Drawn up:**

Vs.	Description of change	Date	Name	Sign
01	First edition	16.06.2010	G. Zimmermann	gez. G. Zimmermann
02	add on risk control measures	05.10.2010	M. Kleiner	gez. M. Kleiner
03	Update to new standard	03.02.2011	G. Zimmermann	<i>G. Zimmermann</i>

Checked and released:

	Division	Date	Name	Sign
Participant 1	R&D	Feb 17, 2011	M. Kleiner	<i>M. Kleiner</i>
Participant 2				
Released by	R&D	Feb 18, 2011	Streit	<i>Streit</i>

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Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D

Identification number: ZI/11/02/03-01

Page 2 of 7

Member code: year/month/day-sequential number

FB 7302 Risk Analysis and control
60601-1-6 MDS SE V5 03 docD

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
5.1. Insufficient labeling of the device (EN 60601-1-6, 5.1, 6.2.2.3) (EN 62366, 4.1, 6)	1. Misapplication	Non or ambiguous labeling	Faulty treatment	Minor	Occasional	ALARP	Labeling of device conforming to standards and well-defined intended use in user's manual	Minor	Inconceivable	Passable	Check list labeling Checklist user's manual
	2. non-ionizing radiation	No or wrong safety precautions	Eye damage Skin damage Fire hazard	Critical	Hardly conceivable	ALARP	Labeling device conforming acc. to standards Warning advice in user's manual	Critical	Hardly conceivable	ALARP	Check list labeling Check list user's manual
5.2. Insufficient or no user's manual (EN 60601-1-6, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong handling and maintenance	Wrong handling Faulty maintenance	Faulty treatment	Minor	Occasional	ALARP	Well-defined description handling and maintenance in user's manual Clear user interface with information for user	Minor	Hardly conceivable	Passable	Check list user's manual
5.1, 5.3 Insufficient labeling and specification of accessories (EN 60601-1-6, 4.2, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. non-ionizing radiation 2. wrong parameters	Uncontrolled emission of non-ionizing radiation Wrong handpiece / spot	Eye damage Overtreatment	Critical	Occasional	ALARP	Labeling acc. to the standards Easy to understand marking of the hand pieces	Critical	Hardly conceivable	ALARP	Check list user's manual Verification of complete device
5.4. Insufficient specification of check-up before using (EN 60601-1-6, 4.2, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong treatment values 2. Uncontrolled emission of non-ionizing radiation 3. Biological contamination	Optic damage Operating in wrong environment Wrong handpiece / spot Use of undesinfected applicator	No therapeutic effect Eye damage Overtreatment Cross infection	Critical	Occasional	ALARP	Description of tests before and during treatment in user's manual Definition of storage-, transportation and operation conditions in user's manual Securing of laser zone by door contact Easy to understand marking of the hand pieces Advice in user's manual	Critical	Hardly conceivable	ALARP	Check list user's manual Verification report 60601-1-4

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision D

Identification number: ZI/11/02/03-01

Member code: year/month/day-sequential number

Page 3 of 7

FB 7302 Risk Analysis and control
60601-1-6 MDS SE vs 03 docD

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
5.5. User's manual too complicated (EN 60601-1-6, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong handling and maintenance	Wrong handling Inadequate maintenance	Faulty treatment	Minor	Occasional	ALARP	Test user's manual by independent persons and users	Minor	Inconceivable	Passable	Check list user's manual Verification complete device
5.6 Insufficient specification for using, operating and maintenance (EN 60601-1-6, 4.2, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong handling	Insufficient description of use	Faulty treatment	minor	Occasional	ALARP	Well-defined description of intended use in user's manual Clear user interface with information for the user	Minor	Hardly conceivable	Passable	Check list user's manual Verification GUI
	2. Deterioration of the equipment qualities	Inadequate maintenance	Faulty treatment Breakdown of the device Harms	Minor	Occasional	ALARP	Well-defined description of maintenance procedure in user's manual Specification of periodically technical safety checks	Minor	Hardly conceivable	Passable	Check list user's manual Verification GUI
5.7. Use by unskilled / untrained personnel (EN 60601-1-6, 4.2.5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Misuse of device	Missing requirements on the user misuse	Insufficient therapeutic effect Faulty treatment	Minor	Occasional	ALARP	Well-defined user profile in user's manual (dentists and medically trained personnel) Securing against unauthorized use by key switch Well-defined description of intended use and handling in user's manual Securing device against unauthorized use by key switch	Minor	Hardly conceivable	Passable	Check list user's manual Verification 60601-1-4
5.8. Well-founded predictable misuse (EN 60601-1-6, 4.2, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. non-intended use	Human error	Injuries	Critical	Occasional	ALARP	Well-defined description of intended use and handling in user's manual.	Critical	Inconceivable	Passable	Check list user's manual Verification 60601-1-4

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D

Identification number: ZI/11/02/03-01

Member code/year/month/day-sequential number

Page 4 of 7

FB 7302 Risk Analysis and control
60601-1-6 MDS SE vs 03 cncD

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
5.9. Insufficient warning of side effects (EN 60601-1-6, 4.2, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong treatment parameters	Not or badly skilled user	Side effects	Minor	Occasional	ALARP	Warning in manuals Instruction of personnel	Minor	Hardly conceivable	Passable	Clinical evaluation Check list Users manual
5.10 Insufficient warning of contraindication (EN 60601-1-6, 4.2.5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	Misuse of the device	Unsufficient description of intended use	Side effects	minor	Occasional	ALARP	Well-defined description of intended use	Minor	Hardly conceivable	Passable	Check list user s manual
6.5. Complex and confusing control and operation system (EN 60601-1-6, 4.2, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Operating error	Too complex control system/ User Interface	Wrong parameter settings Unattended stop of treatment	Minor	Occasional	ALARP	Clear LCD-Display with Touch Screen, self-explaining Laser-Emergency-Stop at the front Blip when change of parameters Ergonomic arrangement of buttons	Minor	Hardly conceivable	Passable	Verification software GUI Check list user s manual
6.7 Arguable or ambiguous display of parameter settings, mea-	1. Unintended or faulty setting of parameters	Confusing display	Non or too minor low therapeutic effect Over treatment Side effects	Critical	Occasional	ALARP	Communication display with confirmation Change of parameters only in Stand-by and check data during change in Ready Self-explaining display	Critical	Hardly conceivable	ALARP	Risk verification 60601-1-4

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D

Identification number: ZI/11/02/03-01

Member code/year/month/day-sequential number

Page 5 of 7

FB 7302 Risk Analysis and control
50671-1-6 MNS SE Vs 03 docD

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
Measurement results or other specifications (EN 60601-1-6, 4.2, 6.1, 6.2.2.3)	2. Misinterpretation of displayed values	Confusing display	Non or too minor therapeutic effect Over treatment, Side effects	Critical	Occasional	ALARP	Clear LCD-Display with Touch Screen Identical positions of the display values in all user-displays	Critical	inconceivable	Passable	Verification software GUI
	3. Confusing display of data	Confusing display	Non or too minor therapeutic effect Over treatment, Side effects	Critical	Occasional	ALARP	Clear LCD-Display with Touch Screen, Self explaining Dispensable data only shown in sub-menus Display of values in SI units	Critical	Inconceivable	Passable	Verification software GUI
6.9. Insufficient visibility, audibility and touch ability (EN 60601-1-6, 4.2, 6.1, 6.2.2.3) (EN 62366 4.1, 5.3)	1. Incorrect setting of parameters	Visibility of display during intensive lightning or through laser goggles	Non or too minor therapeutic effect Over treatment adverse effects	Critical	Occasional	ALARP	Display with backlight laser emergency stop Adjustable brightness of display Clear color design of display, also readable through laser protection glasses	Critical	Inconceivable	Passable	Verification software GUI Verification complete device
unintended break of therapy (EN 60601-1-6, 6.1)	1. unintended release of foot switch	Pressure by foot not high enough	No risk To press foot switch again	Not critical	Occasional	passable	no	NA	NA	No risk	NA
	2. Automatic switch off of the device	Different safety controls of device	No risk	Not critical	Occasional	passable	no	NA	NA	No risk	NA
Amis design (EN 60601-1-6, 6.1) (EN 62366, 4.1, 5.7)	Design of device	Rough edges	hematomas	critical	Occasional	ALARP	design with no rough edges	Minor	Hardly conceivable	passable	Verification complete device

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D

Identification number: ZI/11/02/03-01

Member code/year/month/day-sequential number

Page 6 of 7

FB 7302 Risk Analysis and control
60601-1 & MCS SE V6 03 doc D

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
Repetitive moving (EN 60601-1-6, 6.1) (EN 62366, 4.1, 5.7)	handpiece	Heavy handpiece	Strain of hand	Critical	Occasional	ALARP	weight of handpiece less than 300g	minor	Inconceivable	Passable	Verification complete device
Loud noise (EN 60601-1-6, 6.1) (EN 62366, 4.1, 5.7)	noise	compressor	Damage of hearing	critical	Occasional	ALARP	Compressor with low noise	No risk	Inconceivable	passable	Verification complete device
7.2 No or insufficient maintenance specifications (EN 60601-1-6, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong output values 2. Hazard by electrical shock	Misalignment light or optics Failure electric safety mechanism	No or too minor therapeutic effect electric shock Fire hazard Overtreatment	Critical	Occasional	ALARP	Determination of periodically technical safety check Check temperature Service manual, periodically service training	Critical	Inconceivable	Passable	Check list user's manual Check service manual
7.3 Insufficient maintenance (EN 60601-1-6, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Wrong output values 2. Hazard by electrical shock	Misalignment light or optics Failure electric safety mechanism	Non or to minor therapeutic effect Overtreatment Electric shock Fire hazard	Critical	Occasional	ALARP	Determination of periodically safety check (STK) in users and service manual Instruction personnel when set-up Notice next service date on display Device able to clean easily	Critical	Hardly conceivable	ALARP	Check list user's manual Instruction report Validation complete device

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Revision: D

Identification number: ZI/11/02/03-01

Member code/year/month/day-sequential number

Page 7 of 7

FB 7302 Risk Analysis and control
60601-1-6 MDS SF Vs 03.docD

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
7.4 Missing specification of life time or maintenance intervals (EN 60601-1-6, 5.1, 6.2.2.3) (EN 62366, 4.1, 5.3, 6)	1. Destruction of applicator or optics	Abrasion of applicator tip due to long time use	Non or too minor therapeutic effect Side effects	Minor	Occasional	ALARP	Warning advice in user's manual to change applicator	Minor	Hardly conceivable	Passable	Verification complete device Check list user's manual
8.2 Transport damage (EN 60601-1-6, 4.2, 6.1, 6.2.2.3) (EN 62366, 4.1, 5.3)	1. Too high weight or bad weight spreading 2. Rough edges	Too big, non-ergonomic designed and weighty device, packaging insufficient	1. Injuries during transportation 2. Harms because of weighty moving	Minor	Occasional	ALARP	Weight minimization Evenly distributed weight No squared edges Integrated handle and roles Proper packaging, detectors for shipment	Minor	Inconceivable	Passable	Verification complete device, including packaging

Form 7302	Risk analysis and control acc. to EN ISO 14971:2009	
Revision: D	Identification number: IS/11/01/26-01 Member code/year/month/day-sequential number	Page 1 of 3 FB 7302 Risk Analysis and control 60601-1-2 MDS SE Vs 03.cccD



Medical Product: Mediostar (NeXT) ☐ MDD 93/42/EEC ☐ EN 60601-1-4 1996+ A1:1999 ☐ before R&D
 Item number: 1590 ff ☐ EN 60601-1:2006 ☐ EN 62304:2006 ☒ during R&D
 Assessed construction level: Pre-Serial Device ☐ EN 60601-1-2:2001+ A1:2006 ☐ EN 60601-1-6:2004 ☐ after sales
 Revision user's Manual: 1590 99 001 Vs.05 Jan 11, 2011 ☒ EN 60601-1-2:2007 ☐ EN 60601-1-6:2007 ☐ review reg. changed standards
☐ EN 62366:2008

Related documents:

Drawn up:

Vs.	Description of change	Date	Name	Sign
01	First edition	16.06.2010	M. Kleiner	gez. M. Kleiner
02	add on risk control measures	04.10.2010	M. Kleiner	gez. M. Kleiner
03	Update to new standard	26.01.2011	I. Streit	

Checked and released:

	Division	Date	Name	Sign
Participant 1	PhD	Jan. 26, 2011	Herbst	
Participant 2				
Released by	G. W.	23. Jan. 2011	Thieme	

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Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Asclepion

Laser Technology

Revision: D

Identification number: IS/11/01/26-01

Member code/year/month/day-sequential number


Page 2 of 3

FB 7302 Risk Analysis and control
63601-1-2 MDS SE vs 03 docD

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
1.1 Electricity	1. application part	Induction of energy interference of isolation by electric or magnetic fields	Hazardous voltages	Critical	Inconceivable	Passable	Non	Critical	Inconceivable	Passable	Not applicable
	2. Power supply	electromagnetic effects Harmed isolation	Flash over to signal lines	Critical	Inconceivable	Passable	Non	Critical	Inconceivable	Passable	Not applicable
	3. Source of Supply	Interference of isolation by electric/magnetic fields	Increasing-ground leakage current	Critical	Inconceivable	Passable	Non	Critical	Inconceivable	Passable	Not applicable
1.2 Emitted heat	Internal cooling system	Interference of function by electric or magnetic fields	Overheating	Critical	Inconceivable	Passable	Non	Critical	Inconceivable	Passable	Not applicable
1.5 non-ionizing radiation	1. Internal electronics including electronic shutter	Induction of energy Interference of function by electric or magnetic fields	Unintended change of adjusted parameters	Critical	Hardly conceivable	ALARP	Screening filtering of supply line Electro Static Discharge (ESD)-compatible design	Critical	Inconceivable	Passable	Test report acc. to 60601-1-2
1.13 Magnetic fields	1. Magnetic sensors	Interference of function by external magnetic fields	Generally malfunctions	Critical	Hardly conceivable	ALARP	Use of sensors which are insensitive against external magnetic fields	Critical	Inconceivable	Passable	Check acc. to 60601-1-2

Form 7302

Risk analysis and control acc. to EN ISO 14971:2009

Asclepiion

Revision: D

Identification number: IS/11/01/26-01

Member code/year/month/day-sequential number

Page 3 of 3


FB 7302 Risk Analysis and control
60601-1-2 MDS SE Vx 03 docD



a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
3.1 Electromagnetic noise fields (Immunity)	1 Sensitivity against external EM-noise fields	Insufficient screening against external EMC	Malfunction of the device Uncontrolled emission of non-ionizing radiation pulse / wrong parameters	Critical	Occasional	ALARP	Constructive measurements acc. to EN 60601-1-2 EMC-test of prototype Monitoring of the output parameters	Critical	Inconceivable	Passable	Test report EN 60601-1-2
3.2 Electromagnetic fields (Emission)	1 Disfunction of other devices by EM radiation	Electromagnetic noise fields (grid-bound, air)	Malfunction hazards by malfunction of other devices EMC at user / patients	Critical	Occasional	ALARP	Constructive measurements acc. to EN 60601-1-2, checked at test	Critical	Inconceivable	Passable	Test report EN 60601-1-2
7.1 Faulty data transmission	1 Harmed processor hardware	Induction of energy interference of function by electric or magnetic fields	Wrong output parameters	Critical	Hardly conceivable	ALARP	Screening filtering of supply line Electro Static Discharge (ESD)-compatible design	Critical	Inconceivable	Passable	Check acc. to 60601-1-2
	2 Harmed software	Induction of energy interference of function by electric or magnetic fields	Wrong output parameters	Critical	Hardly conceivable	ALARP	Screening filtering of supply line Electro Static Discharge (ESD)-compatible design	Critical	Inconceivable	Passable	Check acc. to 60601-1-2

Form 7302	Risk analysis and control acc. to DIN EN ISO 14971/ ISO 14971 (in jeweils gültiger Ausgabe entsprechend CL7321 Relevante Sicherheitsnormen)	
Revision: C	Identification number: MMo/11/11/15-01 Member code/year/month/day-sequential number	Page 1 of 8 FB 7302 Risk Analysis and control 93_42 EWG MDS SE V6 04.doc

Medical Product: **Mediostar** ☒ MPR 93/42/EWG ☐ before R&D
 Item number: **1590 - 1599** ☐ DIN EN 60 601-1-2 ☐ during R&D
 Assessed construction level: **Serial Device** ☐ DIN EN 60 601-1-4 ☐ after R&D (before serial production)
 Revision user's Manual: ☐ DIN EN 60 601-1-6 ☒ after production

Drawn up:

Vs.	Description of change	Date	Name	Sign
04	Add on risk control measures (spillage on ME device)	15/11/2011	M.Moller	
03	Add on risk control measures	09/05/2011	M.Moller	gez. M.Moller
02	Add on risk control measures	04/10/2010	M.Kleiner	gez. M.Kleiner
01	First edition	16/06/2010	M.Kleiner	gez. M.Kleiner

Checked and released:	Division	Date	Name	Sign
Participant 1	<i>R&D</i>	<i>Nov 15 2011</i>	<i>Streit</i>	
Participant 2				
Released by	QM	<i>18.10.2011</i>	R.Thieme	

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Revision: C

Identification number: MMo/11/11/15-01

Member code/year/month/day-sequential number

Page 2 of 8

FB 7302 Risk Analysis and control
93 42 EWG MDS SE V8 04 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
1.1. Electricity	1 Power supply	Contact with line voltage accessible parts with mains voltage Failure of basic insulation	Electrical shock Fire	Critical	Hardly conceivable	ALARP	Protective earthing of case with basic insulation to hazardous voltage Double insulation to SIP/SOP Double insulation to applied parts (glass layer) Protective earth of Cooling finger Earthing of casing parts	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
	2 Short circuit in device	Contact with short-circuit current in device and on the case Insulation defect Failure of basic insulation	Fire	Critical	Hardly conceivable	ALARP	Conduction by protective earth conductor 2-pole Main fuse	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
	3.1 Break of protective earth conductor	Contact with hazardous voltage	Electrical shock	Critical	Hardly conceivable	ALARP	Basic insulation Double insulation to SIP/SOP Double insulation to applied parts	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
	3.2 Supply voltage on applied part	External power supply in the case of first fault	Electrical shock	Critical	Hardly conceivable	ALARP	Applied part direct connected with earth (PE) Basic insulation ELV to earth (PE)	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
	4 Moisture on conducting parts	Short circuit in device Contact with line voltage Ingress of moisture (spillage) High humidity (condensation)	Fire	Critical	Occasional	ALARP	Protective insulation Protective earthing Constructive design of case	Critical	Inconceivable	Passable	Test report DIN EN 60601-1 Verification

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Revision: C

Identification number: MMo/11/11/15-01

Member code/year/month/day-sequential number

Page 3 of 8

FB 7302 Risk Analysis and control
93_42 EWG MDS SE V5 04 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
	5. fixing of wiring	Accidental detachment (breaking) of wiring terminal or conductor	Electrical hazards (cable breakage), burning hazards	Critical	Occasional	ALARP	Fixing of wiring and conductors by heat shrink tubes or tie the wires	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
	6. Hazard of damage by edges of wiring / material	dangerous edges	Electrical hazards, burning hazards	Critical	Occasional	ALARP	Construction according to EN 60601-1 Protection of edges (edge protection)	Minor	Inconceivable	Passable	Test report DIN EN 60601-1
	7. Hazard of any means of protection	Wet parts of liquid reservoir (means of protection)	Electrical shock	Critical	Occasional	ALARP	Closed system Filling of water and overflow outside of the device	Critical	Inconceivable	Passable	Test report DIN EN 60601-1 Verification
	8. Hazard of any means of protection	Wet parts of leakage intersections (means of protection)	Electrical hazards	Critical	Occasional	ALARP	Protected earth of housing Leakage (drip) test of all intersections of cooling system	Critical	Inconceivable	Passable	Test report DIN EN 60601-1 Verification
	9. Hazard of any means of protection	Spillage on ME device (means of protection)	Electrical shock	Critical	inconceivable	ALARP	Closed housing system Protected earth of housing No liquid used for application, only US gel in use	Critical	inconceivable	ALARP	Test report DIN EN 60601-1
1.2 Emitted heat	1. Hand piece	Skin burning by touching	Heat up of touchable hand-piece parts by scattered radiation	Minor	Hardly conceivable	ALARP	hand piece test during production T-Sensor in hand piece Warning advice in user's manual	Minor	Inconceivable	Passable	Verification complete device Check list User's manual
	2. Hand piece	Skin burning by touching	Skin cooling damaged	Minor	Hardly conceivable	ALARP	Temperature control: cooling finger cold side, Warning advice in user's manual	Minor	Inconceivable	Passable	verification software Check list User's manual

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Revision C

Identification number: MMo/11/11/15-01

Member code/year/month/day-sequential number

Page 4 of 8


FB 7302 Risk Analysis and Control
53_42 EWG MDS SE V1.04.doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
	3. Wiring/ materials resistance to heat	Burning and melting of material	Burning and insurance	Critical	Occasional	ALARP	Wiring and material parameters according to the standard	Critical	Inconceivable	Passable	CL safety relevant materials, Ball pressure test
1.5. Non-ionizing radiation	1 Non-ionizing radiation	Unauthorized use Incorrect use	Eye damage of patients, user service technician, uninvolved persons Injuries by fright	Critical	Occasional	ALARP	Constructive protection measures in acc. to EN 60825 HA4 (protective case, distant interlock, key-switch, laser emission display) Warning and user labels on device Laser goggles for patients, users and technicians Advice in user's manual Trained staff	Critical	Hardly conceivable	ALARP	verification software CL EN 60825-1 CL EN 60601-2-22 Verification complete device Check-list User's manual
	2 unintentionally release of Non-ionizing radiation	Uncontrolled activation of foot switch	Eye damage Hazards patients, users or uninvolved persons	Critical	Occasional	ALARP	Foot-switch with coverage Standby-Ready Circuit Electronically shutter Wachtdog, Laser-emergency stop Warning advice in user's manual	Critical	Hardly conceivable	ALARP	CL safety relevant materials verification software CL EN 60601-2-22
	3 Laser radiation	Use of laser without eye protection	Eye damage	Critical	Hardly conceivable	ALARP	Laser protection before activating READY Mode Advice in user's manual	Critical	Inconceivable	Passable	Check-list User's manual
	4. Laser radiation after release of foot switch	Footswitch damage, Short circuit signal line to ground potential, Malfunction controller	Eye damage Hazards patients, users or uninvolved persons	Critical	Occasional	ALARP	Foot-switch with two independents switch Electronically shutter, Wachtdog, Laser-Emergency stop	Critical	Hardly conceivable	Passable	Verification complete device

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Revision: C

Identification number: MMo/11/11/15-01

Member code/year/month/day-sequential number

Page 5 of 8

FB 7302 Risk Analysis and control
93_42 EWG MDS SE Vs 04.doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability			Extent of damage	Occurrence probability		
	5. Laser radiation without activation of foot switch	Malfunction controller, Software-Fault	Eye damage Hazards patients, users or uninvolved persons	Critical	Occasional	ALARP	Electronically shutter independent of controller, Hardware control after start Ready-Mode, Laser-Emergency stop, Automatically "Standby-Mode" after inactivity Foot-switch with two independent switch acoustic warning	Critical	Hardly conceivable	Passable	Verification complete device Software Verification
	6. Non-ionizing radiation	Radiation of inflammable or explosive substances Fire or explosion	Combustion smoke poisoning	Critical	Occasional	ALARP	Advice in user's manual Instruction of the staff T-Sensor in handpiece	Critical	Hardly conceivable	ALARP	Check-list User's manual CL safety relevant materials
1.7 Unintended moving	Complete device	Independent action of wheels	Injury and hazard of person and surrounding devices	Minor	Occasional	ALARP	no moveable device, table top	Minor	Inconceivable	passable	Check-list EN 60601-1 User's manual
1.10 Pressure	1. Burst tank (compressor)	Damage of pressure vessels (explosions)	Injury and hazard of person and surrounding devices	Minor	Hardly conceivable	ALARP	Product of pressure and volume is less than 200kPa Internal pressure control (compressor system) Pressure relieve device	Minor	Inconceivable	passable	
	2. active pressure relieve device	Failure of pressure control	Injury and hazard of person and surrounding devices	Minor	Hardly conceivable	ALARP	Metal housing	Minor	Inconceivable	passable	
2.1 Biological contamination	1. Hand-piece	Hand-piece not cleaned	Infection environment caused by contaminated hand-piece	Minor	Hardly conceivable	ALARP	Hand-piece able to disinfect Advice in user's manual Instruction of the staff	Minor	Inconceivable	Passable	Check-list user's manual

Form 7302	Risk analysis and control acc. to DIN EN ISO 14971/ ISO 14971 (in jeweils gültiger Ausgabe entsprechend CL7321 Relevante Sicherheitsnormen)	Asclepion Lasertechnologie
Revision: C	Identification number: MMo/11/11/15-01 Member code/year/month/day-sequential number	Page 6 of 8 FB 7302 Risk Analysis and control 93 42 EWG MDS Se Vs 04 001

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
2.2. Biological incompatibility	1. Infection respectively irritation by contact with applicator	Contact with not biocompatible materials	Infection Irritation	Minor	Hardly conceivable	ALARP	Using of biocompatible materials on contact points	Minor	Inconceivable	Passable	Check-list EN ISO 10993-1
2.5. Allergenicity (hypersensitivity to substances)	1 Hand-piece	Hand-piece not cleaned	Infection of therapy location by contaminated hand-piece	Minor	Hardly conceivable	ALARP	Instruction of the staff Advice in user's manual	Minor	Inconceivable	Passable	Check-list user's manual
2.10 Re- and/or cross infection	1. Hand-piece	Hand-piece not cleaned and sterilized	Infection of therapy location by contaminated hand-piece	Minor	Hardly conceivable	ALARP	Instruction of the staff Advice in user's manual	Minor	Inconceivable	Passable	Check-list user's manual
3.3. Insufficient energy supply	1 Uncontrolled function	Under- or over-voltage	Incorrect therapy caused by incorrect Non-ionizing radiation power	Minor	Occasional	ALARP	voltage-input with rated voltage tolerance of $\pm 10\%$ Monitoring output parameter	Minor	Inconceivable	Passable	Verification complete device Test report DIN EN 60601-1
3.4. Insufficient coolant supply / restriction of cooling	1. Failure of ventilator	Cable break overheating	Incorrect therapy by low Non-ionizing radiation power	Minor	Hardly conceivable	ALARP	Temperature monitoring and shutdown at increased temperatures Constructive secure cable laying	Minor	Inconceivable	Passable	Test report DIN EN 60601-1
	2. Coolant	Blockage, leak in cooling system	Overheating Laserdiode and system	Critical	Occasional	Not passable	- Temperature monitor - Flowsensor - Advice in user's manual	Critical	Hardly conceivable	ALARP	Verification software
3.5. Storage and operating outside the specified ambient conditions	1. Too hot	Too high environmental temperature	Overheating of device Deviation of Non-ionizing radiation parameters	Minor	Occasional	ALARP	Temperature sensor Advice in user's manual	Minor	Hardly conceivable	ALARP	Verification software Check-list user's manual

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Revision: C

Identification number: MMo/11/11/15-01

Member code/year/month/day-sequential number

Page 7 of 8

FB 7302 Risk Analysis and control
93_42 EWG MDS SE V6 04 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
	2. Too moistly	Too high air humidity	Short circuit in device caused by leakage path Fire hazard, Electrical shock	Critical	Occasional	ALARP	Protective earthing Protective insulation of protective conductor	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
	3. Too cold	Ambient temperature under 0°C	Freezing coolant	Minor	Occasional	ALARP	Temperature sensor Advice in user's manual	Minor	Hardly conceivable	ALARP	Verification software CL user manual
	4. too moistly	Humidity outside standard conditions	Electrical breakdown	Critical	Inconceivable	Passable	Only use in house, controlled ambient conditions No risk measures necessary				
3.7. Accidentally mechanic damage	1. Touch of accessible parts, uncontrolled Non-ionizing radiation	Mechanic damage on case	Electrical shock, Eye damage by Non-ionizing radiation, Injuries	Critical	Hardly conceivable	ALARP	Robust case device and handpiece Advice in user's manual Transport package	Critical	Inconceivable	Passable	Check-list user's manual Test report EN 60601-1
5.15. Arras and square edges	Hazard of injuries by edges, cant over	Defects in construction, production	Injuries of mechanics or patients	Minor	Hardly conceivable	ALARP	Construction according to EN 60601-1	Minor	Inconceivable	Passable	Test report DIN EN 60601-1
5.3 Insufficient specs of two applied accessories	Exchange of wrong accessories	Use not suitable accessories	Electrical shock, Eye damage by Non-ionizing radiation Injuries	Minor	Hardly conceivable	ALARP	Special hand piece interface (mechanical and electrical) Code protected hand piece instructions in user's manual	Minor	Inconceivable	Passable	Verification of laser Checklist user manual
6.6. Arguable or in-explicit status of device	1. Uncontrolled function	Selfactivation or malfunction of user interface, control system	Incorrect therapy fire	Critical	Hardly conceivable	ALARP	Electrical and EMC-safe construction of device Check function when switch on and switch ready Laser-EMERGENCY-STOP, key switch and all-pole main switch for hazardous shut down	Critical	Inconceivable	Passable	Test report DIN EN 60601-1 Test report DIN EN 60601-1-2

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Revision: C

Identification number: MMo/11/11/15-01

Member code/year/month/day-sequential number

Page 8 of 8

 FB 7302 Risk Analysis and control
 93 42 EWG MGS SE Vs 04 doc

a) Hazards	b) Hazard source	c) Potential causes	d) Potential consequences / risks	e) Risk assessment before the precaution		f) Risk assessment	g) Technical safety precautions to be taken	h) Risk assessment after the precaution		i) Residual risk assessment	j) Proof
				Extent of damage	Occurrence probability	Passable ALARP Not passable		Extent of damage	Occurrence probability	Passable ALARP Not passable	
7.5. Loss electrical/mechanical integrity	1. Hazard by electrical shock or mechanical exposure or emitted Non-ionizing radiation	Transport damage Failure of safety equipment	Electric shock Injuries Eye damage by laser radiation	Critical	Hardly conceivable	ALARP	Construction according to EN60601-1	Critical	Inconceivable	Passable	Test report DIN EN 60601-1
7.8. Deterioration of function as result of repeated using	1. Incorrect output values 2. electric shock	Wrong adjustment laser or optics Failure of safety equipment	Electrical shock non or too minor therapeutic effect	Critical	Occasional	ALARP	Determination of annual safety technical control in the user's manual and service manual instruction of the staff when setting up	Critical	Hardly conceivable	ALARP	Check-list user's manual users briefing protocol
7.9. Breakdown device / function by failure / breakdown component part	Failure of automatic shutoff uncontrolled Non-ionizing radiation emission	Failure of laser protection device Failure of safety device Failure transmission system	Electric shock eye damage by laser radiation	Critical	Occasional	ALARP	electronic shutter Watch-dog system Determination of annual safety technical control in the user's manual and service manual instruction of the staff when setting up Regular control of homogenous beam pattern by burn paper	Critical	Inconceivable	Passable	Verification software Check-list user's manual Users briefing protocol
8. Breakdown air clearance	Failure of air clearance insulation	Too less barometric pressure	Electrical shock	Critical	Occasional	ALARP	Designed to fulfill the requirements of patient protection till an altitude of 3000m (standard requirements of 3 rd edition 60601-1)	Critical	Inconceivable	Passable	Check by inspection

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices		 Seite 1 von 15 CL 73523 EN 62366- MDS 159x Vs 01 ZI-2011-03-31 01.rtf
	Anderungs- Index: A	Identifikationsnummer: ZI/2011/03/31-01 Kürzel/Jahr/Monat/Tag lfd. Nr.	

Date of issue: 31.3.11

Total number of pages.....: 15

Testing Laboratory: Asclepion Laser Technologies GmbH

Address: Brüsseler Straße 10
07747 Jena

Tested by:

29.3.11 Zimmermann / B. Zimmermann
Date / Name / Signature

Released by:

31.03.11 Kaiser, D. Kaiser
Date / Name/ Signature

Test specification:

Standard: IEC 62366:2007(ed.1)/ EN 62366: 2008

Test procedure.....: CB

Non-standard test
method.....: N/A

Test Report Form No.: IEC62366A

Test Report Form Originator.....: TÜV Rheinland Product Safety GmbH

Master TRF: Dated 2008-07

This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.

Test item description.....: Laser device

Trade Mark: MeDioStar SE 159x
(advancement of MeDioStar XT/miXT 153x
Proof of usability of the existing device:
FB 73511 Review 2011 MeDioStar XT-miXT ZI-2011-03-28-02)

Manufacturer: Asclepion Laser Technologies

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Testing procedure and testing location:	
☐ CB Testing Laboratory:	
Testing location/ address.....:	
☐ Associated CB Test Laboratory:	
Testing location/ address.....:	
Tested by (name + signature)..: Approved by (+ signature).....:	
☐ Testing procedure: TMP	
Tested by (name + signature)..: Approved by (+ signature).....:	
Testing location/ address.....:	
☐ Testing procedure: WMT	
Tested by (name + signature)..: Witnessed by (+ signature).....:	
Approved by (+ signature).....:	
Testing location/ address.....:	
☐ Testing procedure: SMT	
Tested by (name + signature)..: Approved by (+ signature).....:	
Supervised by (+ signature):	
Testing location/ address.....:	
☐ Testing procedure: RMT	
Tested by (name + signature)..: Approved by (+ signature).....:	
Supervised by (+ signature):	
Testing location/ address.....:	

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: Kürzel/Jahr/Monat/Tag-Idd. Nr	Seite 3 von 15 CL 73523 EN 62366 AID: 159x Vs.01 ZI-2011 05-01 01 df

Summary of testing:

Tests performed (name of test and test clause):	Testing location: Asclepion Laser Technologies Brüsseler Str. 10 D-07747 Jena Germany
---	---

Summary of compliance with National Differences: Conformity evaluation according to EN 62366:2008

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: Kurz/Jahr/Monat/Tag-Idd. Nr	Seite 4 von 15 CL 73523 EN 62366:2008 159x Vs 01_ZI-2011_0431- 01.rtf

Copy of marking plate

Asclepion Laser Technologies GmbH

 07747 Jena, Brüsseler Straße 10, Germany

Model: 159X

MedioStar

	1	2	3	4	5	6	7	8	9	0
X										X

S/N: 159X00102

 12/2010

1275

100-240 Vac

50 / 60 Hz

1500 VA

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 5 von 15 CL 73523 EN 62366 MD 159x Vs 01 ZI-2011-03-31 01.nf

Test item particulars.....:

Classification of installation and
use.....:

Supply connection
.....:

Context of Use
.....:

Abbreviations used in the report:

- Usability Engineering:	UE	- Risk analysis:	RA
- User interface:	UI	- Risk management:	RM
- Primary operating function:	POF		

Possible test case verdicts:

- test case does not apply to the test object : N/A
- test object does meet the requirement : Pass (P)
- test object does not meet the requirement : Fail (F)

Testing:

Date of receipt of test items
.....:

Date(s) of performance of tests :

General remarks:

The test results presented in this report relate only to the object tested.
This report shall not be reproduced, except in full, without the written approval of the Issuing testing laboratory.

"(see Enclosure #)" refers to additional information appended to the report.
"(see appended table)" refers to a table appended to the report.

Throughout this report, a point (coma) is used as the decimal separator.

List of test equipment must be kept on file and available for review.

This Test Report contains the general safety requirements as related to the usability of Medical Equipment.

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Änderungs- Index A	Identifikationsnummer.: Kürzel/Jahr/Monat/Tag lfd. Nr	Seite 6 von 15 CL 73523 EN 62366 MOS 199x Vs 01 Zi 2011-03-31- 01.rtf

General product information:

Specification	MeDioStar
Device model	Floor-based unit
Display	8.4" LCD display
Operator guidance	Touchscreen
Cooling	Internal cooling cycle
Door interlock connector	5 V / 10 mA
Laser warning lamp	Potential-free relay contact max. 24 V / 1 A (make contact)
Permissible ambient conditions	Temperature: 15 °C to 30 °C rel. air humidity: max. 70 % (no condensation) height: max. 2000m above sea level
Permissible transport and storage conditions	Temperature: 0°C to +50°C rel. air humidity: 10 % to 95 % (no condensation)
Dimensions	365 x 560 x 465 mm (W x D x H)
Weight	without handles and handpiece holder ca. 35 kg
Power requirements	100 -240VAC, 50/60 Hz, max. 1500 VA
Main fuse	overcurrent release 16A, medium slow-blow
Overvoltage category	II (IEC 60664-1)
Max. peak of mains voltage	V MT peak = 2500V peak
Max. nominal mains voltage	V MN r.m.s. = 300V r.m.s.
Classification according to RL 93/42/EEC	II b
Operating mode	The device has been designed for continuous operation

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Änderungs-Index: A	Identifikationsnummer.: Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 7 von 15 CL 73523 EN 62366 MDS 159x Vs 01 ZI-2011 03 31- 01.rtf

Specification

MeDioStar

Laser model	High-power diode laser array
Wavelength	800 - 950 nm
Laser class	4
Required laser safety goggles (acc. To DIN EN 207: 2009)	800 - 950 nm D LB5
Pulse duration	max. 400 ms (depending on the mode and handpiece)
Pulse energy	max. 60 J
Pulse frequency	max. 12 Hz (depending on the mode and handpiece)
Spot sizes (handpieces)	standard: 13x10 mm ² (STD) optional: 4x3 mm ² (VAS) 8x6 mm ² (HP)
Tolerance of output power	± 20 %
Laser beam mode	multimode
Laser beam diameter at the handpiece output	see spot sizes
Energy density	max. 44 J/cm ² @ spot size 13 x 10 mm ²
	max. 90 J/cm ² @ spot size 8 x 6 mm ²
	max. 210 J/cm ² @ spot size 4 x 3 mm ²
Beam divergence at the handpiece output (round angle, 1/e ²)	typ. 400 mrad @ handpiece STD typ. 700 mrad @ handpiece HP typ. 600 mrad @ handpiece VAS
Nominal ocular hazard distance (with handpiece)	35 m
Protection class of the part used	B

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 8 von 15 CL 73523 EN 62366 MDS 159x Vs.01 ZI 2011/03/31 01.rtf

4	GENERAL REQUIREMENTS		
4.1	General Requirements		
4.1.1	USABILITY ENGINEERING PROCESS		
	Has the MANUFACTURER established, documented and maintained a USABILITY ENGINEERING PROCESS to provide SAFETY for the PATIENT, USER and others related to USABILITY for the product?	FB 7302 Risc analysis and control 60601-1-6 MDS SE Vs 03 – ZI/11/02/03-01-Konformitätsakte MeDioStar 002/10, 7.2.2	P
	Does the PROCESS addressed USER INTERACTIONS with the MEDICAL DEVICE according to the ACCOMPANYING DOCUMENT including, but not limited to transport, storage, installation, operation, maintenance, repair and disposal?	User Manual 1590 99 001 (deutsch) 1590 99 002 (english) Konformitätsakte MeDioStar 002/10, 2.4	P
4.1.2	RESIDUAL RISK		
	Are RESIDUAL RISKS associated with USABILITY of the MEDICAL DEVICE presumed to be acceptable, unless there is OBJECTIVE EVIDENCE to the contrary and documented?	FB 7302 Risc analysis and control 60601-1-6 MDS SE Vs 03 – ZI/11/02/03-01-Konformitätsakte MeDioStar 002/10, 7.2.2	P
4.1.3	Information for SAFETY		
	MANUFACTURER subject the information for safety used as a RISK CONTROL to the USABILITY ENGINEERING PROCESS (e.g., warnings or limitation of use in the ACCOMPANYING DOCUMENTS, marking, etc.) in the UE file and Accompanying documents.	User Manual 1590 99 001 (deutsch) 1590 99 002 (english) Konformitätsakte MeDioStar 002/10, 2.4 incl. labeling as specified in this document	P
	Disregarding such information for SAFETY is considered beyond any further reasonable means of RISK CONTROL	FB 7302 Risc analysis and control 60601-1-6 MDS SE Vs 03 – ZI/11/02/03-01-Konformitätsakte MeDioStar 002/10, 7.2.2	P
4.2	USABILITY ENGINEERING FILE		
	The results of the USABILITY ENGINEERING PROCESS are recorded in the USABILITY ENGINEERING FILE	FB 7351 ValidiERG Gebrauchstgkt MDS NeXT Vs.01 ZI/2011/03/29-1 Konformitätsakte MeDioStar 002/10, 7.6	P

CL 73523

IEC 62366/ EN 62366 : 2008

Medical devices

Application of usability engineering to medical devices

Anderungs-
Index: A

Identifikationsnummer.: ZI/2011/03/31-01

Kurz-/Jahr/Monat/Tag-lfd. Nr.

Seite 9 von 15
CL 73523 EN 62366: MED
159x Vs 01 ZI-2011-3-3
01.nf

	The records and other documents that make up the USABILITY ENGINEERING FILE form part of other documents and files (e.g., a MANUFACTURER'S product design file or RISK MANAGEMENT FILE), (SEE List of documents make up the UE file)	Usability Engineering file is a part of the Konformitätsakte MeDioStar 002/10 – Chapter 7.6
4.3	Scaling of the USABILITY ENGINEERING effort	
	The USABILITY ENGINEERING PROCESS is scaled based on the significance of any modifications depending on the results of the RISK ANALYSIS and documented	FB 7302 Risc analysis and control 60601-1-6 MDS SE Vs 03 - ZI/11/02/03-01- Konformitätsakte MeDioStar 002/10, 7.2.2
5	USABILITY ENGINEERING PROCESS	
5.1	Application specification	
	Application of MEDICAL DEVICE in the USABILITY ENGINEERING FILE is specified by the MANUFACTURER and includes	Application not changed MeDioStar SE 159x is an advancement of MeDioStar XT/miXT 153x Proof of usability of the existing device: FB 73511 Review 2011 MeDioStar XT-miXT ZI-2011-03-28-02
	– intended medical indication (e.g., conditions(s) or disease(s) to be screened, monitored, treated, diagnosed, or prevented);	User Manual 1590 99 001 (deutsch) 1590 99 002 (english) Konformitätsakte MeDioStar 002/10, 2.4
	– intended PATIENT population (e.g., age, weight, health, condition);	
	– intended part of the body or type of tissue applied to or interacted with;	
	– intended conditions of use (e.g., environment including hygienic requirements, frequency of use, location, mobility); and	
	– operating principle(s)	
5.2	Frequently used functions	
	Are frequently used functions that involve USER interaction with the MEDICAL DEVICE are determined and recorded in the USABILITY ENGINEERING FILE?	Like MeDioStar 1530 according to requirement specification FB 7321 Lastenheft MeDioStar SE ZI/10/11/18-1

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 10 von 15 CL 73523 EN 62366:2008 159x Vs 01 ZI/2011/03/31 01.rtf

5.3	Identification of HAZARDS and HAZARDOUS SITUATIONS related to USABILITY		
5.3.1	Identification of characteristics to SAFETY		
	Identification of characteristics related to SAFETY (part of a RISK ANALYSIS) that focuses on USABILITY performed according to ISO 14971:2007, 4.2.	FB 7302 Risc analysis and control 60601-1-6 MDS SE Vs 03 – ZI/11/02/03-01-Konformitätsakte MeDioStar 002/10, 7.2.2	P
	During the identification characteristics related to SAFETY, the following are considered: – application specification, including USER PROFILE(S); and – frequently used functions.		P
	Results of this identification characteristics related to SAFETY recorded in the USABILITY ENGINEERING FILE		P
5.3.2	Identification of known or foreseeable HAZARDS and HAZARDOUS SITUATIONS		
	MANUFACTURER has identified known or foreseeable HAZARDS (part of a RISK ANALYSIS) related to USABILITY according to ISO 14971:2007, 4.3.		P
	Identification of HAZARDS considered HAZARDS to PATIENTS, USERS and other persons	FB 7302 Risc analysis and control 60601-1-6 MDS SE Vs 03 – ZI/11/02/03-01-Konformitätsakte MeDioStar 002/10, 7.2.2	P
	Reasonably foreseeable sequences or combinations of events involving the USER INTERFACE that can result in a HAZARDOUS SITUATION associated with the MEDICAL DEVICE were identified. The SEVERITY of the resulting possible HARM is determined.		P

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Änderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 11 von 15 CL 73523 EN 62366 MDC 159x Vs.01 ZI/2011/03-31 01.rtf

	During the identification of HAZARDS and HAZARDOUS SITUATIONS, the following was considered: – application specification, including USER PROFILE(S); – task related requirements; – context of use; – information on HAZARDS and HAZARDOUS SITUATIONS known for existing USER INTERFACES of MEDICAL DEVICES of a similar type, if available; – preliminary USE SCENARIOS; – possible USE ERRORS; – if an incorrect mental model of the operation of the MEDICAL DEVICE can cause a USE ERROR resulting in a HAZARDOUS SITUATION; and – results of the review of the USER INTERFACE		P
	The results of this identification of HAZARDS, HAZARDOUS SITUATIONS and SEVERITY are recorded in the USABILITY ENGINEERING FILE.		P
5.4	PRIMARY OPERATING FUNCTIONS		
	The manufacturer has determined the PRIMARY OPERATING FUNCTIONS and recorded in the USABILITY ENGINEERING FILE	Appendix 1a of FB 7351 Validierg Gebrauchstgkt MDS NeXT Vs.01 ZI-2011-03-29-1 Konformitätsakte MeDioStar 002/10, 7.6	P
	The inputs to the PRIMARY OPERATING FUNCTIONS include frequently used functions and functions related to SAFETY of the MEDICAL DEVICE		P
5.5	USABILITY SPECIFICATION		
	MANUFACTURER developed a USABILITY SPECIFICATION recorded in the USABILITY ENGINEERING FILE as part of the USABILITY ENGINEERING PROCESS	according to requirement specification FB 7321 Lastenheft MeDioStar SE ZI/10/11/18-1 Konformitätsakte MeDioStar 002/10,5.1	P
	The USABILITY SPECIFICATION recorded in USABILITY ENGINEERING FILE. The USABILITY SPECIFICATION may be integrated into other specifications		P

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 12 von 15 CL 73523 EN 62366:2008 159x Vs.01 ZI 2011 02/01 01.rst

	The USABILITY SPECIFICATION includes: – application specification; – PRIMARY OPERATING FUNCTIONS – HAZARDS and HAZARDOUS SITUATIONS related to the USABILITY; and – known or foreseeable USE ERRORS associated with the MEDICAL DEVICE		P
	The USABILITY SPECIFICATION describes at least: – USE SCENARIOS related to the PRIMARY OPERATING FUNCTIONS EQUIPMENT, including – frequent Use Scenarios, and – reasonably foreseeable worst case USE SCENARIOS;		P
	– USER INTERFACE requirements for the PRIMARY OPERATING FUNCTIONS EQUIPMENT including those to mitigate Risk;		P
	– Requirements for determining whether PRIMARY OPERATING FUNCTIONS are easily recognizable by the USER.		P
5.6	USABILITY VALIDATION plan		
	The MANUFACTURER has developed and maintains a USABILITY VALIDATION plan specifying:	The validation was focused to the changes to the MeDioStar XT/miXT 153x FB 7351 ValidiERG Gebrauchstgkt MDS NeXT Vs.01 ZI/2011/03/29-1 Konformitätsakte MeDioStar 002/10, 7.6	P
	– any method used for VALIDATION of the USABILITY of the PRIMARY OPERATING FUNCTIONS;		P
	– the criteria for determining successful VALIDATION of the USABILITY of the PRIMARY OPERATING FUNCTIONS based on the USABILITY SPECIFICATION; and		P
	– the involvement of representative intended USERS		N/A
	USABILITY VALIDATION performed in a laboratory setting		N/A

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Anderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kurz-/Jahr/Monat/Tag-lfd. Nr	Seite 13 von 15 CL 73523 EN 62366 MDS 159x Vs.01 ZI/2011/03/31 01/11

	USABILITY VALIDATION performed in a simulated use environment.....:		P
	USABILITY VALIDATION performed in the actual use environment.....:	--	N/A
	The USABILITY VALIDATION plan addresses: – frequent Use Scenarios, and – reasonably foreseeable worst case USE SCENARIOS that are identified in the USABILITY SPECIFICATION	The validation was focused to the changes	P
	The USABILITY VALIDATION plan recorded in the USABILITY ENGINEERING FILE		P
5.7	USER INTERFACE design and implementation		
	MANUFACTURER designed and implemented the USER INTERFACE as described in the USABILITY SPECIFICATION utilizing, as appropriate, USABILITY ENGINEERING methods and techniques	Proof of usability of the existing device: FB 73511 Review 2011 MeDioStar XT-miXT ZI-2011-03-28-02 FB 7351 ValidiERG Gebrauchstgkt MDS NeXT Vs.01 ZI/2011/03/29-1 Konformitätsakte MeDioStar 002/10, 7.6	P
5.8	USABILITY VERIFICATION		
	MANUFACTURER verified the implementation of the MEDICAL DEVICE USER INTERFACE design according to the USABILITY SPECIFICATION	FB 7351 Verifizierung MDS SE Vs.01 110308 MMO/11/03/08-01 Konformitätsakte MeDioStar 002/10, 7.9	P
	The results of the verification are recorded in USABILITY ENGINEERING FILE		P
5.9	USABILITY VALIDATION		
	The MANUFACTURER has validated the USABILITY of the MEDICAL DEVICE according to the USABILITY VALIDATION plan	FB 7351 ValidiERG Gebrauchstgkt MDS NeXT Vs.01 ZI/2011/03/29-1 Konformitätsakte MeDioStar 002/10, 7.6	P
	The results are recorded in the USABILITY ENGINEERING FILE		P

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	Seite 14 von 15 CL 73523 EN 62366 IUDS 159x Vs. 01 ZI-2011-03-31 01 rtf
Anderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kurzel/Jahr/Monat/Tag-Ild. Nr	


	For the acceptance criteria documented in the USABILITY VALIDATION plan that are not met: - further USER INTERFACE design and implementation activities are performed; or - if further improvement is not practicable, the MANUFACTURER may gather and review data and literature to determine if the medical benefits of the INTENDED USE outweigh the RISK arising from USABILITY problems To perform this step, the MANUFACTURER needs to estimate the RISK arising from USABILITY problems.		N/A
6	ACCOMPANYING DOCUMENTS		
	The ACCOMPANYING DOCUMENT includes a summary of the MEDICAL DEVICE application specification		P
	A concise description of the MEDICAL DEVICE, its operating principles, significant physical and performance characteristics and intended USER PROFILE are included in the ACCOMPANYING DOCUMENT		P
	The ACCOMPANYING DOCUMENT is written at a level consistent with the intended OPERATOR PROFILE	User Manual 1590 99 001 (deutsch) 1590 99 002 (english) Konformitätsakte MeDioStar 002/10, 2.4	P
	The ACCOMPANYING DOCUMENT for equipment are, optionally, provided electronically		N/A
	USABILITY ENGINEERING PROCESS includes the information that will need to be provided as a hard copy or as markings on MEDICAL DEVICE when ACCOMPANYING DOCUMENTS are provided electronically		N/A
7	TRAINING AND MATERIALS FOR TRAINING		
	The required training on the MEDICAL DEVICE for safe and effective use of PRIMARY OPERATING FUNCTIONS by the intended USER is given by:		N/A
	- necessary training materials provided by the manufacturer;	No special training necessary	N/A
	- necessary training materials are available; or		N/A
	- the manufacturer provides TRAINING		N/A

CL 73523	IEC 62366/ EN 62366 : 2008 Medical devices Application of usability engineering to medical devices	
Änderungs- Index: A	Identifikationsnummer.: ZI/2011/03/31-01 Kürzel/Jahr/Monat/Tag-lfd. Nr	Seite 15 von 15 CL 73523 EN 62366 MDS 159x Vs 01 ZI 2011-03-31 01 df

	The ACCOMPANYING DOCUMENT describes the available training options (Recommendation: ACCOMPANYING DOCUMENT include the suggested duration and frequency of such training)		N/A
	INTENDED USE AND USER PROFILE(S) are the basis for TRAINING and TRAINING material		N/A

List of documents make up the UE file

4.2	Document title	Version
	CL 73523 EN 62366 MDS 159x Vs.01 ZI-2011-03-31-01 ZI/11/03/31-01	01
	CL 7364 EN 1041 MDS 159x Vs.01 IS/11/03/15-02	01
	CL 7364 EN 60601-1-6 MDS 159x Vs.01 ZI-2011-03-30-01	01
	FB 7351 Validlerg Gebrauchstgkt MDS NeXT Vs.01 ZI-2011-03-29-1	01
	Anl 1a zu FB 7351- Val Gebrtgkt MDS - Hauptbedien POFs ZI-07-03-12-1	01
	Anl 1b zu FB 7351 - Val Gebrtgkt MDS - Software-Ergonomie-Akte	01
	Anl 1c zu FB 7351 - Val Gebrtgkt MDS - Zusammenfassung Evaluierung GUI	01
	Anl 2 zu FB 7351 - Val Gebrtgkt MDS - Parameterfeld 03-11	01
	FB 73511 Review 2011 MeDioStar XT-miXT ZI-2011-03-28-02	01
	FB 73511 Review 2011 Anl.3 ZI-11-03-28-01	01
	FB 7321 Lastenheft MeDioStar SE Vs02 ZI/10/11/18-1	02
supplementary information: Files for risc analysis and control see chapter 7.2.2 of Konformitätsakte MeDioStar 002/10, 7.2.2		

CL 7360	Conformity Evaluation acc. to EN 60 601-1-6: 2007	 Asclepion Laser Technologies
Revision index: C	Identification number: ZI/2011/03/30-01 Identification code/year/month/day-number	page 1 of 14 C: 7360 EN 60601-1-6 MDS 159x Vs.01 ZI-2011-03-30-01.doc

Tested equipment: MeDioStar
Name and/or model: 159x x = 1...9
Manufacturer: Asclepion Laser Technologies GmbH
Address: Brüsseler Str. 10, 07747 Jena, Germany


Applied standard: EN 60 601-1-6: 2007 (DIN EN 60 601-1-6: 2008- 02)
 Medical electrical equipment – general requirements for safety - collateral standard: Usability

Attachments:

Examined by: Gabi Zimmermann
 Name Date

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Asclepion
Test TechnologyRevision index:
C

Identification number: ZI/2011/03/30-01

Identification code/year/month/day-number

page 2 of 14
CL 7360 EN 60601-1-6 MOS 159x
Vs 01 ZI 2011-03-30-01 doc

Identification

Details of examiner:

Identification of the person responsible for examining and classifying the product under inspection.

Name: Gabi Zimmermann
print full name (block letters)Position: Senior Product Manager
print full title

If the above named person is not an employee of the manufacturer of the product to be examined, state the details of the examiner's employer or organization:

Organisation: _____

Address: _____


Medical product

Details of product to be inspected:

User's manual 1590 99 001_1-28. Feb. 2011

☒ portable equipment ☐ mobile equipment ☐ permanently installed equipment

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Revision index:
C

Identification number: ZI/2011/03/30-01

Identification code/year/month/day number

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Revision index:
C

Identification number: ZI/2011/03/30-01

Identification code/year/month/day-number.

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



Examined by: M. Möller June 18, 2014
 Name Date Sign 
Checked by:  June 18, 2014
 Name Date Sign 

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Identification <i>Details of examiner:</i> Identification of the person responsible for examining the product under inspection: Name: Michael Moller print full name (block letters) Position: Test engineer R&D print full title If the above named person is not an employee of the manufacturer of the product to be examined, state the details of the examiner's employer or organization. Organisation: _____ Address: _____ _____ Medizinprodukt/Medical device Details of product to be inspected: User's manual: 159099001-Vs04-5.Febr2014-GA-MSPRO-DE	<i>Picture of the device</i> 
---	---

CL 7351	Grundlegende Anforderungen zur Bewertung der Konformität nach RL 93/42/EWG * Conformity Evaluation Essential Requirements acc. to CD 93/42/EEC *	 page 3 of 33 CL 7351 93_42_EEC MSNEXT V06 MM0140518-01 docx
Revision index: C	Identification number: MMo/14/06/18-01 Identification code/year/month/day-number	

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RL 93/42/EWG *

Conformity Evaluation Essential Requirements acc. to CD 93/42/EEC *

ASCIEM

Revision index: C

Identification number: MMo/14/06/18-01

Identification code/year/month/day-number

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FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month-sequential no.	Page 1 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01.doc

Verification	☐	Validation	☒
Product / part:	MeDioStar SE	Material number:	159x
Cause:	☐ Product development	☒ Product change	☐
Original document:	☐ Requirement specification of Lastenheft ZI/10/11/18-01	☐ Design specification of	☐
	☐ Risk analysis of	☐ Change-application no.	☐
Serial number:	N/A		

Brief description: (traceability of testing): literature route

The goal of the new version of the MeDioStar was to improve the technique, which is some years old and to reduce the production costs. Therefore there was no change relating to the intended use as described in the User Manual of the MeDioStar so far. This was the request of the Requirement Specification (Lastenheft ZI/10/11/18-01) and was realized in the development of the new version. There is no need for a clinical study.

This assessment includes the evaluation of actual literature of lasers with the same specs like the laser MeDioStar and of our own experience from our customers of the MeDioStar versions.

Used measuring and test equipment: not applicable

Documents of the enclosure: (test reports etc.)



Requirement Specification (Lastenheft ZI/10/11/18-01), Konformitätsakte MeDioStar ALT 002/10, 5.1
 App. 1: Protocol for Review of Literature
 App. 2: Literature: Clinical Studies with comparable devices and own experience with MedioStar so far (including assessment of parameters used in the studies)
 App. 3: App. 3 to Reviews 2011
 App. 4: Overview diode lasers used for clinical studies

Final assessment of results: (brief assessment, summary measures)

Based on the realized risk analysis and the evaluation of the information arising by the analysis of the literature the final evaluation, that the benefit of the new device and the application is higher than the risk for the patient and the user. The expected side effects and risks are marginal and not higher than the risks of other devices (including the older versions of the MeDioStar) which are used since years for adequate applications.

☐ VALIDATION NEGATIV ☒ VALIDATION POSITIV

if negative, derived measures

Plan drawn by:			Validation performed:		
6.1.11	G. Zimmermann		6.10.11	G. Zimmermann	
Date	Name	Signature	Date	Name	Signature
Plan checked and released:			Validation checked and released:		
7.1.11	A. Steps		7.10.11	A. Steps	
Date	Name	Signature	Date	Name	Signature
Cause of inconsistency (if first issue) Drawn by: G. Zimmermann R&D: 1.07.05 J. Stroh 			FB checked and released by QA: 10.09.05 R. Thiele 		

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no	Page 2 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01.doc

Clinical Evaluation MeDioStar SE

Contents

- 0 Preliminary note
- 1 Qualification of the author and the reviewer
- 2 Argument for literature route
- 3 Evaluation of literature
 - 3.1 Study design
 - 3.2 Evaluation of the scientific studies regarding effectiveness
 - 3.2.1 Effectiveness demonstrated in percent
 - 3.2.2 Effectiveness demonstrated by a scale
 - 3.2.3 Comparison of results
 - 3.2.4 histological findings
 - 3.3 Evaluation of the scientific studies regarding side effects
 - 3.3.1 Evaluation of side effects reported in the studies in detail
 - 3.3.2 Side effects evaluated by histological findings (a) and objective measurements (b)
 - 3.4 Evaluation of contraindications/exclusion criteria described in the studies
 - 3.4.1 Exclusion criteria of the studies and reports
 - 3.4.2 Contraindications of the consulted standard book
 - 3.5 Evaluation of customer experience
 - 3.6 Evaluation of reports of the BfArM
- 4 Final assessment of results

0 Preliminary Note

The high power diode lasers in the near infrared (NIR) wavelength between 800 and 950nm have now being used for more than 10 years with very good results. Compared to other lasers and intense pulse light systems (IPL) for the same indications the diode lasers have a lower tendency to side effects.

**Table of laser devices and light sources
which are used for the intended use of the MeDioStar :**

Type of device	Wavelength	comment
Ruby lasers	694nm	Relatively high risk of side effects specially for darker skin because of high melanin absorption
Alexandrite lasers	755nm	Relatively high risk of side effects specially for darker skin because of high melanin absorption
High power diode lasers	800 – 950nm	Low risk of side effects because of lower melanin absorption
Nd:YAG lasers (long pulsed)	1064nm	Low risk of side effects because of low melanin absorption. More pain during treatment compared to all other devices because of high water absorption. Water is the main component of the skin.
IPL	Range between 600 – 1200nm dependent on device	Highest risk of side effects because of broad range of wavelengths including a part of visible light, which has the highest melanin absorption compared to the lasers

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month/ sequential no.	Page 3 of 13 FB 7351e ClinEval MS-SE ZI-2011 10 05-01.doc

Asclepion Laser Technologies developed the first version of MeDioStar in 1999 (name of the company to that time was Aesculap-Meditec) and the production and sales started in 2000. Since that time up to now 1476 high power diode lasers were sold (in different MeDioStar versions) worldwide, including countries where patients live with mainly dark skin.

The customers treated their patients successfully and with very low side effects. These facts were already shown in the literature of the clinical evaluation of the last MeDioStar version. The overview of these articles is attached as Part A (studies, which were already evaluated for the clinical evaluation of the current version of MeDioStar) of Appendix 2 of this evaluation.

Among the other studies, specially study no. 2, which is a comparison between MeDioStar to the only other high power diode laser to that time 'LightSheer' evaluated², that with same efficacy the MeDioStar has less erythema and burned hair (no pigment changes or scars for both). Additionally there was a customer experience with MeDioStar (no. 4) with a very high number of patients / treatments (200 / 978). The statistical analysis of these treatments of patients with all skin types showed a percentage of side effects of 0.5% only (transient hyperpigmentation).

The following Clinical Evaluation of the MeDioStar SE is therefore based on the experience with the current version of MeDioStar and analyzes additionally the literature which was not yet included in the clinical evaluation of the current MeDioStar version.

1 Qualification of the author and the reviewer

The author of this evaluation, Mrs. G. Zimmermann is the Senior Product Manager and the person in duty for medical validation and application by Asclepion Laser Technologies with a lot of experience to carry out clinical evaluations.

Mrs. G. Zimmermann has a diploma in physics and has experience in application of medical lasers for more than 20 years, especially for surgical lasers, lasers for dermatology and aesthetics including the MeDioStar and devices using intense pulsed light (IPL), radiofrequency and acoustic waves.

After graduation she developed laser devices as a scientific engineer for more than 10 years. Since 1988 she was responsible for medical application of different lasers including MeDioStar and worked together with physicians in various specialties. She also performed laser safety courses including medical application for physicians for 18 years now and educated the sales staff and doctors in using the devices.

She was in charge to supervising own application studies and evaluating scientific literature. Working as Product Manager for more than 20 years she was also responsible for fixing of application parameters of new lasers as well as supporting the development in respect to ergonomically questions, writing the application part of user manuals, evaluating new devices in respect to medical use and ergonomics.

The reviewer of this evaluation, Mr. André Steps, is not an employee of the manufacturer Asclepion Laser Technologies GmbH.

Mr. Steps has a prediploma in medicine and a sport science degree. Since more than 2 years he educates doctors and nurses in the theoretical principles and practical application of lasers, IPLs and other devices for dermatological and esthetic treatments. In this time he did a vast number of trainings all over the world.

Therefore he has an excellent overview about all relevant lasers and devices used worldwide for the intended use of the MeDioStar as well

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	 Asclepion Laser Technologies
Revision Index A	Identification number: ZI/2011/10/05-01 Member Code: year/month/sequential no	Page 4 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01.doc

2 Argument for literature route

The goal of the new version of the MeDioStar was to improve the technique, which is some years old, and to reduce the production costs.

The intended use (including indications and target group/contraindications) of the new version MeDioStar SE is identical compared to the current version of the MeDioStar. The requested application parameters related to the clinical performance, benefit/risk ratio, safety and risks are the same. These are the range of wavelength, of fluence and of pulse duration, which were used, as well as the skin cooling system.

The reason for this requirement in the Requirement Specification (Lastenheft ZI/10/11/18-01) was, that the effectivity of treatments with the MeDioStar was proven, also compared to similar devices. Equally important, during all the time there were no complaints about serious damage to health. The number of side effects (small blisters, crusts and temporary pigment changes), which cannot be excluded totally, was very low.

Additionally the remaining risks are identical to the current MeDioStar, which did not cause any harm to patients or users in the past 12 years besides in some cases the possible side effects stated in the manual. Therefore the warnings in the User Manual are the same. All warnings and informations necessary in the User Manual of the new device are checked for completeness at the basis of the Risk analysis MK 10/10/04-01 and MK 10/10/05-01 by CL7356 Checkliste Gebrauchsanweisung, ZI/11/03/04-01 (see Record of Conformity MeDioStar 1590 ALT 002/10 chapter 7.2 and 7.9).

As for the current version of MeDioStar the treatment of the indications of the intended use by high power diode lasers (including parameters) are standard treatments. There were already scientific studies for the clinical evaluation of the current MeDioStar version showing the efficacy and safety of these treatments for the patients.

Of course in the meantime more studies are published with the MeDioStar and similar lasers, which support the positive evaluation some years ago.

That's why it is possible to demonstrate, that the use of the new MeDioStar SE for the intended use in question will be effective and safe with a reasonable benefit/risk ratio for the patients treated.

Evaluating these facts there is no need for a new clinical study. The assessment by literature route is possible.

This assessment includes the evaluation of actual literature of lasers with the same specs like the laser MeDioStar and of our own experience from our customers of the MeDioStar versions.

3 Evaluation of literature

This assessment includes the evaluation of actual literature of lasers with the same specs like the laser MeDioStar and of our own experience from our customers of the MeDioStar versions, which are sold in the past 10 years.

The following appendixes are part of this clinical evaluation:

App. 1: Protocol for Review of Literature

App. 2: Literature.

Clinical Studies with comparable devices and own experience with MeDioStar so far
 (including assessment of lasers/parameters used in the studies)
 (including contraindications stated in a standard book for intended use in question)

App. 3: App. 3 to Reviews 2011 Usability (including clinical evaluation)

App. 4: Overview diode lasers used for clinical studies

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index A	Identification number: ZI/2011/10/05-01 Member Code/year/month/ sequential no.	Page 5 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01 doc

3.1 Study design

- (A) All studies (no. 10 – 19 of App. 2) confirm to the inclusion criterions as described in App. 1 and are prospective studies.

Note:

The reports no. 20, 21, 22 and 23 are not considered as a 'study'. But it was included in this evaluation because it shows additionally the good results of our customers (see 3.5 for details).

- (B) All studies demonstrate the effectiveness of the treatments in question as well as the safety with a reasonable benefit/risk ratio
- (C) All studies were published in journals with very good reputation, despite no. 10.

The study no. 10 was not published in a journal and is already some years old (because it was not published, it was not known so far). The reason to include this paper in this evaluation was because of 2 reasons (additionally to the conformity to the inclusion criterions as for all studies):

1. The co-author Rox Anderson, MD
He has the best reputation for lasers used for dermatological and esthetic treatments worldwide. He published a vast number of scientific articles and even standard books for lasers in dermatology including the scientific basic principles for these applications. For about 20 years now the community follows his principles. Additionally he educated a lot of doctors, who are now well-known experts in this field.
 2. The laser used in this study (brand: LightSheer/ company now: Lumenis)
This device is the most similar device to the MeDioStar at the market and is sold (in different versions, but with the same intended use and parameters) as long as the MeDioStar.
This laser is the 'gold standard' for the laser epilation, because it was sold with the highest numbers all over the world and was proven for effectiveness and safety.
Already in the clinical evaluation of the current MeDioStar version it was shown by study (App. 2: Part A, no. 2) that the MeDioStar has the same effect and less side effects compared to the LightSheer.
- (D) All clinical studies were realized with the MeDioStar (no. 15 and 16 of App. 2) or similar devices (no. 10, 11, 12, 13, 14, 17, 18, 19). The comparability of the other lasers used in the studies is demonstrated in the assessments of the appendix 2. The range of parameters used in the studies can be realized with the MeDioStar as well.
- (E) Assessment of the study design in detail:
(see also the column Study Design of App. 2)

Study 10

See also comment under (C) above
valuable design, because:

- different combinations of application parameters were tested at test sites and compared under the same conditions of 40 patients
- the same parameter combinations are already used, but so far not directly compared in this way
- long term follow up (up to 6 months), which is necessary for hair removal
- all skin types were included.

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no	Page 6 of 13 FB 7351e ClnEval MS-SE ZI-2011-10-05-01 doc

This is not easy to do, because dependent on location normally not all skin types are available. That's why very seldom all skin types are included in a study. (This is important, because different skin types have a different tendency to side effects.)

Study 11

Very valuable design, because

- comparison between the LightSheer (see comment C) and other lasers used for the same indications
- the other lasers are not exactly comparable to the MeDioStar because of other wavelengths, but are competition at the market
- High number of patients (171)
- Very long follow up (12 months)

Study 12

Very valuable design, because

- direct comparison between 2 high power diodes like MeDioStar (LightSheer and the other brand Soprano/ company Alma Lasers)
- Direct comparison between 2 different techniques, which are both used for hair removal.
highest tolerable fluence with 1 pass vs. low fluence with multiple passes
- Direct comparison between light skin (skin type I-II) and dark skin (III-IV)
- Assessment for different hair (color and thickness)
- Long follow up with 6 months

Study 13

Valuable design, because

- direct comparison between a high power diode like MeDioStar (brand: Soprano/ company: Alma Lasers) vs. an Alexandrite laser, which is also a very well recognized laser for epilation
- disadvantage: 18 patients only (but several areas)
- Long follow up with 6 months

Study 14

Very valuable design, because

- Multicenter study
- Very high number of patients (368)
- Objective measurement of biophysical properties of the skin additionally to effectiveness and side effects
- Long follow up with 6 months

Study 15

Very valuable design, because

- Very scientific, randomized study with MeDioStar
- Comparison: treated vs. untreated (control) axilla
- Objective measurement of biophysical properties of the skin additionally to effectiveness and side effects

Study 16

Very valuable design, because

- study with MeDioStar
- high number of patients (143)

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no	Page 7 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01.doc

- treatment of Asian skin (Chinese)
(important, because reaction of Asian skin is different compared to all other skin types)

Study 17

Valuable design, because

- randomized, blinded split-face study
- Comparison: 940nm diode vs. 532nm laser, which is also well accepted in the market for vascular treatments (brand: VariLite / company: Iridex)

Study 18

Valuable design, because

- Histological assessment of 2 biopsies
- Disadvantage: 15 patients (but multiple vessels)

Study 19

Valuable design, because

- Vascular as well as pigmented lesions in the face were treated
- The additional use of a chromophor for a part of the lesions plays no role for this evaluation

3.2 Evaluation of the scientific studies regarding effectiveness

All these new studies document as well as the older studies of the last MeDioStar evaluation the efficacy of treatments with the MeDioStar or similar devices for the intended use.

It is already known, that the result is individually different and that is also shown in these studies, e.g. the result of hair reduction depends par example on hair type. This is therefore explained in the User Manual of the MeDioStar (chapter Application).

For details see chapters below:

3.2.1 Effectiveness demonstrated in percent:

The reduction of hair for the diode laser used was shown in percent as follows:

- 6 months after a single treatment with high fluence: 28 – 50%
higher fluences have better effectiveness (no. 10)
- 6 and 12 months after multiple treatments: 65% (no. 11)
- 6 months after multiple treatments (no. 12):
86% with low fluence, multiple passes
91% with high fluence, one pass
- 6 months after 3 treatments (no. 13): 2 – 43 %
(the 2% are for blond and thin hair, which is known to have bad results)
- 6 weeks after a single treatment: 48% (no. 15)

Improvement of ectatic vessels by the diode laser in percent was shown as:

- 2 months after 2nd treatment: between 57,1% (chin) and 70% (cheeks) (no. 17)
- 2 months after 2nd treatment: 55.1% (no. 18)
- 4 weeks after 1 treatment all vascular lesions had healed (no. 19)

Improvement of skin relating to pigments (rejuvenation) was shown as:

- 2 months after 2nd treatment all pigmented lesions had healed (no. 19)

3.2.2 Effectiveness demonstrated by a scale:

In study no. 14 the degree of clearance of hair density was shown as follows:

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no.	Page 8 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01.doc

reduction of hair between	0 – 25%:	29 patients
	25 – 49%:	102 patients
	50 – 74%:	219 patients
	75 – 100%:	18 patients

3.2.3. Comparison of results

- for hair reduction the diode had the best results compared to Alexandrite and Ruby laser, which had the poorest results (no. 11)
- result for hair reduction of one pass of diode laser with highest tolerable fluence was slightly better compared to low fluence of diode laser with multiple passes (see 3.2.1 also, no. 10)
- for dark skin types the diode laser with low fluence was clearly better compared to Alexandrite laser, whereas the Alexandrite was slightly better for light skin (no. 13)
- for vascular treatments the diode laser was significantly more efficacious for larger caliber vessels compared to the 532nm laser, both wavelengths equally for smaller vessels

3.2.4 histological findings:

- cytopathic and vacuole changes, also focal epidermolysis and architectural changes in the hair shaft were shown (no. 14), which is relevant for hair removal
- immediately after treatment of veins the histological assessment showed that the veins had collapsed, completely surrounded by mildly coagulated tissue (no. 18)

3.3 Evaluation of the scientific studies regarding side effects

All these new studies document as well as the older studies of the last MeDioStar evaluation the very low risk of unwanted side effects of treatments with the MeDioStar or similar devices for the intended use with a reasonable benefit/risk ratio for the patients. 2 new studies which include histological findings (no. 9) and objective measurements of biophysical skin properties (no. 15) document, that these treatments are safe.

All side effects seen with the 838 patients treated in all studies (much more sessions!) were transient, besides 2 slight hypochromic halos in case of 2 re-treated pigments. Pigment changes, if appeared, often need more than 2 months to disappear. Side effects are in general more likely the darker the skin type of the patient is (additional Asian skin) and for pigment removal and vascular treatments because of higher absorption in these cases.

The size of the side effects in case of pigment removal (as reported above) or vessel treatment is normally very small because limited to the size of the lesion treated. Side effects cannot be totally excluded, but reduced to a very low risk by doing test treatments. This is explained in the User Manual of the MeDioStar (chapter Application). For details see below:

3.3.1. Evaluation of side effects reported in the studies in detail:

1. Transient erythema and perifollicular edema were mentioned in study no. 14 for most patients and sessions.

Comment by author of the evaluation: This is seen as the desired clinical end point for the test treatment – therefore not seen as real side effect in other studies and that's why often not mentioned

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no	Page 9 of 13 FB 7351e ClinEval MS SE ZI-2011 10-05 01 doc

2. In all studies no permanent side effects were observed with the diode lasers (such as scars or permanent hypo-/hyperpigmentation, which could be possible in general) despite of a slight hypochromic halo in case of the 2 re-treated melanoses after 2 months (reported in no. 19 for pigment removal).

Especially the statistics of the study no. 14 is very valuable, because out of 1840 treatments of patients with dark skin (higher risk of side effects!) realized by 3 centers only transient side effects were seen (5.5%). This was mainly pseudofolliculitis, which is often seen even after normal shaving in dark skinned people.

3. No side effects at all were reported for the diode laser in study no. 13 (epilation, 18 patients, 54 treatment sessions) and in study no. 17 (vascular treatments, 24 facial sites treated, 48 sessions)

4. In comparison to other lasers the least complications were observed with the diode laser:

- compared to Alexandrite and Ruby (no. 11, epilation, 171 patients)
- compared to Alexandrite (no. 13, epilation, 18 patients, 54 sessions)
no pain, no side effects with diode -
1 case of crusts and pigment change as well with Alexandrite
- compared to 532nm (no. 17, vascular, 24 facial sites, 48 sessions)
less pain, erythema, swelling, crusting with 532nm only

5. Temporary pigment changes using diode lasers were reported in study no. 10 (see below), few cases in no. 11 (lowest number for diode compared to other lasers) and 7 cases out of 143 patients (4.9%) with Asian skin (higher risk of side effects!)

All patients of the darkest skin type VI had transient pigment change in study no. 10. The fluence used of 15J/cm² was because of the study design and normally too high for this skin type. Their recommendation at the end of the study was, to start the test treatment for skin type VI with a fluence of 10 J/cm² - that's reasonable.

6. One minor burns of a patients leg was reported (no. 12 and no. 18)

7. Less pain with the diode laser was reported in no. 13 and no. 17

3.3.2. Side effects evaluated by histological findings (a) and objective measurements (b):

- a) Epidermis, stratum corneum and keratin were intact with no identifiable changes immediately post treatment (no. 14, epilation, 9 biopsies).
- b) No changes in Melanin index, trans-epidermal water loss, capacity of stratum corneum hydration, sebum and pH level up to 6 weeks after treatment with diode laser (no. 15, epilation, 31 patients).

3.4 Evaluation of contraindications/exclusion criteria described in the studies

In the evaluation of the contraindications one has to consider, that not all exclusion criteria of the studies are real contraindications. Some of them are included only e.g. because the result of the study has to be clearly related to the laser treatment in question (point 11) or because of the special study design (e.g. related to special skin types, one have to exclude the others, e.g. no 14).

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no	Page 10 of 13 FB 7351e ClinEval MS-SE ZI-2011 10-05-01.doc

That's why additionally to the studies (listed in 3.1) and reports (listed in 3.4) a standard book, which describes application of different lasers in dermatology including contraindications was consulted (details see 3.3.2)

The exclusion criteria are in detail listed and evaluated in 3.3.1 for the studies mentioned, the contraindications of the standard book in 3.3.2.

All together one can clearly see, that the contraindications of the experts are of course not as large-scaled compared to the list of exclusion criteria of the studies.

Nevertheless all real exclusion criteria of the studies and contraindications of the standard book are included in the user manual of the MeDioStar.

3.4.1 Exclusion criteria of the studies and reports

The points 1 to 4 and 6 to 9 are stated as contraindications in the user manual of the MeDioStar (chapters 1.3, 6.4.4 and 6.5.2). Skin disease in the treatment area of point 8 is for Herpes also stated in the chapters mentioned before. Open wounds, infection and inflammation (point 8) are sometimes related to spider leg veins and are therefore stated in the manual as contraindications particularly for the treatment of spider leg veins, specified as ulcer, vasculitis and thrombophlebitis.

The exclusion criteria 'Nursing' (no. 5) is not seen as a contraindication, because there is no additional possible side effect compared to other treatments (additionally: not mentioned as contraindication in the standard book, see 3.4.2).

The exclusion criteria 'hormonal dysfunction' (no. 10) is related to the known effect, that a high value of androgenetic hormones in women stimulates more hair growth. This would of course influence the outcome of a study, that's why these patients were excluded in study no. 12. The fact, that the result of the hair removal treatment (number of treatment sessions) depends among others on the hormonal status is specified in chapter 6.4.10 of the User Manual.

List of exclusion criteria of the different studies:

1. Photosensitivity and medication with photosensitivity as known side effect:
stated in studies no. 12, 15 and 22
2. Tanned skin
stated as exclusion criteria in study no. 12, 13, 15, 18 and 22
3. Cancerous and pre-cancerous lesions
stated in studies no. 12, 19 and 22
4. Pregnancy
stated as exclusion criteria in studies no. 12, 14, 15, 17 and 22
5. Nursing
stated as exclusion criteria in studies no. 12, 14, 15 and 18
6. Children
stated as exclusion criteria in studies no. 14 and 15
7. Keloids or hypertrophic scar formation
stated as exclusion criteria in studies no. 12, 13, 14, 15, 17, 19 and 22
8. Skin disease or inflammation, open wounds, infection, Herpes
stated as exclusion criteria in studies no. 12, 13, 14, 15, 17, 18 and 22

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no	Page 11 of 13 FB 7351e ClinEval MS-SE ZI 2011-10-05 01.doc

9. Tattoos in the treatment area stated as exclusion criteria in studies no. 12 and 22
10. Hormonal dysfunction or intake of medications with androgenetic effects stated as exclusion criteria in study no. 12
11. Other exclusion criteria such as e.g. other epilation methods used before, which are clearly stated to meet the study design (no. 12 – 15, 17)

3.4.2 Contraindications of the consulted standard book

Consulted standard book (see also App. 2, Part B):

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Both authors are experts for application of lasers in dermatology and esthetic medicine in Germany educating newcomers.

In this book among others the contraindications are listed for the Alexandrite (Kardorff) and the Nd:YAG laser (Drosner) for the intended use in question.

For the diode laser one has to consider, that the wavelength of the diode laser MeDioStar (800-950nm), is between Alexandrite (755nm) and Nd:YAG (1064nm). The melanin absorption, which is important for the possible side effects with respect to pigment changes and therefore for the contraindications, is lower with longer wavelength. As one can see in the list below for treatments with the Alexandrite laser (1) more caution is necessary (B, C, D) compared to Nd:YAG (2). This fact was shown also for the diode laser similar to the MeDioStar e.g. in comparison study no. 11.

Therefore it is not necessary to observe C and D for the MeDioStar. The points A, B, E, F, G (cancerous and pre-cancerous lesions in the treatment area) and I are included in the contraindications of the MeDioStar user manual.

Point H is not really a danger (not mentioned for Alexandrite or in all other studies).

The content of point J is described in the chapter 6.4.3 Selection of patients of the user manual.

For vascular treatments point K is considered in the user manual by the indication limited to vessel size between 0.5 and 2mm (chapter 6.5.1 Selection of patients) and point L as contraindication 'thrombophlebitis'.

List of contraindications

- A. Intake of medication which causes photosensitivity (1)
- B. Tanned skin (1)
- C. Tendency to hypo- or hyperpigmentation (1)
- D. Skin type V and VI (1)
- E. abnormal scar formation (1)
- F. anamnesis of Herpes (prophylaxis) (1)
- G. nevi (2)
- H. with deliberation treatment of large areas of the head (2)
- I. strict indication for treatment of children (2)
- J. gray (1 + 2) or blond (2) hair and fine hair (2)

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no.	Page 12 of 13 FB 7351e ClinEval MS-SE ZI-2011-10-05-01.doc

additionally for vascular treatments:

- K. very fine vessels (redness of cheeks and in the face) (1)
- L. coagulopathy (2)

3.5 Evaluation of customer experience

Scientific studies realized with MeDioStar (see above no. 15 and no. 16) and reports of customer experience (see below no. 20 – 23) demonstrate the effectiveness and safety with a reasonable benefit/ risk ratio using the current version of MeDioStar for the intended use.

Additionally the direct comparison of the wavelength of the MeDioStar versions (no. 22 and no. 23) demonstrate once more clearly that the difference relating to efficacy is not statistically significant and the side effects are very low for both with a slight advantage for the 940nm for dark skin types, which was the reason to add this wavelength into the MeDioStar some years ago.

For details see below:

Report 20

This report was included, because a high number of patients (106) were treated specially for vascular treatments and a long follow-up was realized (6 months).

The effectiveness was reported as follows:

- complete healing – 75% (79 patients),
- significant improvement – 20% (21 patients),
- moderate improvement – 4 % (4 patients),
- inconsiderable improvement – 1% (1 patient).

Report 21

This report demonstrates by an objective measurement with a special camera the improvement of the skin after rejuvenation treatment. Despite only 6 patients were treated once, this report was included, because the improvement of pore size, wrinkles and texture is usually done by evaluation of patient satisfaction and of before and after photos by an observer or the therapist. That is not objective.

The camera with a special software (Canfield Visia System) assesses the improvement in a percentage scale by measurements of before and after photos.

The results after only 1 treatment were (3 treatments are recommended):

- Improvement of skin texture: 22 – 31%
- Improvement of wrinkles: 33 – 39%
- Improvement of pore size: 10 – 19%
- Improvement of age spots: 15%

Report 22

In this report 5 patients were treated using a side by side comparison. It was included despite of the low number of patients, because this is a direct comparison with the same patient, same device from all technical details and identical parameters besides the wavelength.

The result for these patients with mainly light skin types was:

- no notable differences between the wavelengths in efficiency and discomfort
- no side effects were seen.

FB 7351e	Verification / Validation Report Clinical Evaluation MeDioStar SE	
Revision Index: A	Identification number: ZI/2011/10/05-01 Member Code/year/month/sequential no.	Page 13 of 13 FB 7351e: ClinEval MeDioStar SE 10-05-21 doc

Report 23

As with report no. 22 the same design was used for the treatments, but in this case with 25 patients including dark skin types.

The result was:

- no significant difference in effectiveness between the wavelengths
average hair reduction for all patients: 58.11% for 810 and 60.35% for 940nm
- only 1 patients had 20% difference between both sides
- no side effects with 940nm, 2 side effects with 810nm
(1 crust, 1 prolonged erythema and hyperpigmentation)
- specially the dark skinned patients felt the 940nm as being more comfortable

3.6 Evaluation of reports of the BfArM (German Federal Department for Drugs and Medical Devices)

The last research of the recommendations of the BfArM as well as the recalls of manufacturers was done at the 28th of March 2011 and gave no results which would be relevant for the MeDioStar.

For details see App. 3 of this evaluation.

4 Final assessment of results

This evaluation of studies and reports realized with the older version of MeDioStar or other devices which were demonstrated as equivalent showed the efficacy of the treatments indicated in the intended use of this new MeDioStar laser device.

Side effects due to individual overtreatment can be blisters, crusts, pigment changes and scars, which are small, because limited to the target (e.g. hair, vessel). They are normally transient despite the small scars, which are seldom. Side effects cannot be totally avoided, but as shown in the evaluation the diode lasers and the MeDioStar cause less side effects compared to other lasers. Additionally the side effects can be reduced to a minimum by doing a test treatment as described in the user manual.

The contraindications found in the literature are included in the manual.

The user is informed about remaining risks in the user manual:

- the possible side effects in chapters 6.4.5 and 6.5.3
- the contraindications in chapters 1.3, 6.4.4 and 6.5.2


Based on this analysis of the literature the final evaluation, that the benefit of the new device and the application is higher than the risk for the patient and the user. The expected side effects and risks are marginal and not higher than the risks of other devices (including the older versions of the MeDioStar) which are used since years for the applications in question.

G. Zimmermann

I hereby certify that above copy is a true copy of the original document which has been produced to me.

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