CryoPro

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CryoPro

User Manual

Cryosurgical instrument CryoPro in the following manufacturing modifications:

CryoPro Maxi, CryoPro Mini, with accessories

Manufactured by Cortex Technology ApS, Denmark

Copyright©2015 Cortex Technology ApS

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Contents

1. Intended use of the cryosurgical instrument CryoPro	3
2. General characteristics of the cryosurgical instrument	3
2.1. Description and pictures of the cryosurgical instrument	3
2.2. General parameters of the cryosurgical instrument and its accessories	7
3. Indications for use of the medical device CryoPro	17
4. Counter-indications	17
5. Safety measures	18
6. Pre-starting and operating procedures of the CryoPro device	20
6.1. General directions	20
6.2. Lesion freezing procedure	20
6.3. Closed contact probes	21
6.4. Cervical probes	26
6.5. Open spray tips	28
6.6. Additional tips (adapters)	31
6.7. Protective accessories	34
7. Storage precautions for the cryosurgical instrument	36
8. Maintenance, disinfection and disposal of the cryosurgical instrument	36
9. Transportation of the cryosurgical instrument	37
10. Complete delivery set of the cryosurgical instrument	37
10.1. Cryosurgical instrument CryoPro, manufacturing modification CryoPro Maxi	37
10.2. Cryosurgical instrument CryoPro, manufacturing modification CryoPro Mini	38
11. Packing and marking of the cryosurgical instrument	39
12. Operational lifetime of the cryosurgical instrument	43
13. Warranty	43
13.1. Conditions of use of warranty	43
14. Requirements on technical servicing and repair of the medical device CryoPro ...	43
15. Manufacturer of the medical device CryoPro	43
16. Authorized representative of the manufacturer	44

1. Intended use of the cryosurgical instrument CryoPro

Cryosurgical instrument CryoPro in the following manufacturing modifications: CryoPro Maxi, CryoPro Mini, with accessories is intended to be used for removing lesions on skin and mucous membranes by their exposure to the temperature of the liquid nitrogen.

2. General characteristics of the cryosurgical instrument


2.1. Description and pictures of the cryosurgical instrument

Cryosurgical instrument CryoPro is produced in the following manufacturing modifications:


Cryosurgical instrument CryoPro Maxi
(Picture 1, Page 4),

Cryosurgical instrument CryoPro Mini
(Picture 2, Page 5).


The basic set of the cryosurgical instruments CryoPro Maxi and CryoPro Mini includes a set of open spray tips (Pictures 3 and 4, Page 6).


Picture 1. Cryosurgical instrument CryoPro Maxi.

Picture 2. Cryosurgical instrument CryoPro Mini.

Picture 3. A set of open spray tips and a stand.

Picture 4. A set of open spray tips.


- | | |
|---------------------------|---------|
| 1. Open spray tip - OS BS | - 1 ea |
| 2. Open spray tip - OS A | - 1 ea |
| 3. Open spray tip - OS B | - 1 ea |
| 4. Open spray tip - OS C | - 2 ea |
| 5. Open spray tip - OS D | - 1 ea |
| 6. Stand for spray tips | - 1 ea. |

2.2. General parameters of the cryosurgical instrument and its accessories

Table 1. Operational parameters of the cryosurgical instrument

	Manufacturing modification	
Characteristics	CryoPro Maxi	CryoPro Mini
Useful capacity of the can	0,5 l	0,35 l
Dimensions of the cryosurgical instrument	Height 287,0±1,0 mm, Diameter of the can 68,0±0,5 mm	Height 223,0±1,0 mm, Diameter of the can 68,0±0,5 mm
Weight of the cryosurgical instrument (without cooling agent)	Not more than 660 g	Not more than 610 g
Weight of the cryosurgical instrument (with cooling agent)	Not more than 1100 g	Not more than 900 g
Operating temperature of the tips	Down to -120°C	Down to -120°C
The period of static storage of liquid nitrogen in the cryosurgical instrument	20 to 24 hours	12 to 14 hours
Cooling agent used	Liquid nitrogen at t°= -196°C	
Working pressure in the can	Not more than 0,7 bar	
Liquid nitrogen yield rate of open spray tips: • type OS A • type OS B • type OS C • type OS D • type OS LL • type OS SS • type OS BS • type OS SOFT	34,8 g/min 19,5 g/min 11,0 g/min 4,9 g/min 19,5 g/min 8,7 g/min 8,7 g/min 23,7 g/min	

Table 2.1. Size and weight of accessories: Closed Contact Probes


- | | |
|--------------------------------------|--|
| 1 – Length of the probe | 6 – Length of the working part of the head |
| 2 – Diameter of the knurled body | 7 – Diameter of the working part of the head |
| 3 – Length of the tube of the body | 8 – Length of the waste pipe |
| 4 – Diameter of the tube of the body | 9 – Outer diameter of the waste pipe |
| 5 – Length of the tip's head | 10 – Inner diameter of the waste pipe |

Description	Length of probe, mm	Diameter of the knurled body, mm	Length/ diameter of the tube of the body, mm	Length of the head/ length/ diameter of the working part of the probe's head, mm	Length/ outer diameter/ inner diameter of the waste pipe, mm	Weight, g
Closed contact probe - CP SP	67,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	9,5±0,1/ 5,2±0,1/ 6,5±0,1	250,0 /5,0/3,0	14,0±1,0
Closed contact probe – CP 1 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,0±0,1/ 3,0±0,1/ 1,0±0,1	250,0 /5,0/3,0	13,0±1,0
Closed contact probe – CP 2 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,5±0,1/ 3,2±0,1/ 2,0±0,1	250,0 /5,0/3,0	13,0±1,0
Closed contact probe – CP 3 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,5±0,1/ 3,5±0,1/ 3,0±0,1	250,0 /5,0/3,0	13,0±1,0
Closed contact probe – CP 4 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,8±0,1/ 3,7±0,1/ 4,0±0,1	250,0 /5,0/3,0	14,0±1,0

Closed contact probe – CP 5 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,8±0,1/ 4,2±0,1/ 5,0±0,1	250,0 /5,0/3,0	14,0±1,0
Closed contact probe – CP 6 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	9,2±0,1/ 5,0±0,1/ 6,0±0,1	250,0 /5,0/3,0	14,0±1,0
Closed contact probe – CP 8 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	7,5±0,1/ 2,0±0,1/ 8,0±0,1	250,0 /5,0/3,0	15,0±1,0
Closed contact probe – CP 10 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	7,5±0,1/ 2,2±0,1/ 10,0±0,1	250,0 /5,0/3,0	16,0±1,0
Closed contact probe – CP 15 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,0±0,1/ 2,7±0,1/ 15,0±0,1	250,0 /5,0/3,0	17,0±1,0
Closed contact probe – CP 20 mm	65,0±0,2	10,0±0,1	41,0±0,1/ 6,0±0,1	8,0±0,1/ 2,8±0,1/ 20,0±0,1	250,0 /5,0/3,0	20,0±1,0

Table 2.2. Size and weight of accessories: Cervical Probes


- 1 – Length of probe
- 2 – Diameter of knurled body
- 3 – Length of the tube of the body
- 4 – Diameter of the tube of the body
- 5 – Length of the probe's head
- 6 – Diameter of the probe's head

- 7 – Length of the head's finger
- 8 – Diameter of the heads' finger
- 9 – Length of the waste pipe
- 10 – Outer diameter of the waste pipe
- 11 – Inner diameter of the waste pipe

Наименование	Length of probe, mm	Diameter of the knurled body, mm	Length/ diameter of the tube of the body, mm	Length/ diameter of the probe's head, mm	Length/ diameter of the probe's finger, mm	Length/ outer diameter/ inner diameter of the waste pipe, mm	Weight, g
Cervical probe 19x5 mm	186,0±0,2	10,0±0,1	156±1,0/ 6,5±0,5	17,0/ 19,0	5,0/5,0	460,0/5,0/3,0	31,0±1,0
Cervical probe 19x15 mm	197,0±0,2	10,0±0,1	156±1,0/ 6,5±0,5	30,0/ 19,0	15,0/5,0	460,0/5,0/3,0	33,0±1,0
Cervical probe 25x5 mm	189,0±0,2	10,0±0,1	156±1,0/ 6,5±0,5	19,0/ 25,0	5,0/5,0	460,0/5,0/3,0	39,0±1,0
Cervical probe 19 mm, dome shaped	184,0±0,2	10,0±0,1	156±1,0/ 6,5±0,5	12,0/19,0	--	460,0/5,0/3,0	32,0±1,0
Cervical probe 25 mm, dome shaped	185,0±0,2	10,0±0,1	156±1,0/ 6,5±0,5	14,0/25,0	--	460,0/5,0/3,0	39,0±1,0

Table 2.3. Size and weight of accessories: Straight and Bent Open Spray Tips


1 – Length of the tip
2 – Diameter of the knurled body

3 – Length of the tube of the body
4 – Diameter of the tube of the body
5 – Diameter of the outlet of the tube

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Table 2.4. Size and weight of accessories: Open Spray Tips OS A, OS B, OS C, OS D

Open Spray Tip OS A.


1 – Length of the tip

3 – Diameter of the orifice

2 – Diameter of the knurled body

Description	Length of the tip, mm	Diameter of the knurled body, mm	Diameter of the orifice, mm	Weight, g
Open spray tip - OS A	11,0±0,2	10,0±0,1	0,8±0,05	3,6±0,5
Open spray tip - OS B	11,0±0,2	10,0±0,1	0,6±0,05	3,6±0,5
Open spray tip - OS C	11,0±0,2	10,0±0,1	0,45±0,05	3,6±0,5
Open spray tip - OS D	11,0±0,2	10,0±0,1	0,3±0,05	3,6±0,5

Table 2.5. Size and weight of accessories: Open Spray Tip - OS SOFT


1 – Length of the tip
2 – Diameter of the knurled body

3 – Outlets

Description	Length of the tip, mm	Diameter of the knurled body, mm	Number / diameter of the outlets in mm	Weight, g
Open spray tip - OS SOFT	25,0±0,2	10,0±0,1	9 orifices $\varnothing 0,22 \pm 0,05$	9,5±1,0

Table 2.6. Size and weight of accessories: Open Spray Tip - OS LL

A – Open Spray Tip OS LL, disassembled view

B – Open Spray Tip OS LL, work-ready view.


1 – Length of the tip
2 – Diameter of the tip
3 – Length of the brass tip body
4 – Length of the body of the stainless steel Luer Lock adapter

5 – Diameter of the outlet of the body of the Luer Lock adapter
6 – Length/diameter/pitch/type of the inner threading of the brass tip's body
7 – Length/diameter/pitch/type of the fitting of the Luer Lock adapter

Description	Length/ diameter of the tip, mm	Length of the brass tip's body, mm	Length of the body of the Luer Lock adapter, mm	Diameter of the outlet of the body of the Luer Lock adapter, mm	Length/ diameter/ pitch/ type of the inner threading of the brass tip's body, mm	Length/ diameter/ pitch/ type of the fitting of the Luer Lock adapter, mm	Length/ diameter/ type of knurling of the body of the Luer Lock adapter, mm	Weight, g
-------------	---------------------------------	------------------------------------	---	---	--	---	---	-----------

Open spray tip - OS LL	28,0±0,1 / 10,0±0,1	18,5±0,1	18,5±0,1	2,3±0,05	9,0±0,05/ 6,0/ 0,8±0,05/ M6	9,0±0,05/ 6,0/ 0,8±0,05/ M6	22±0,05/ 10,0±0,05/ reticular/ pitch 1,5±0,1	15,0± 1,0
------------------------	---------------------	----------	----------	----------	-----------------------------	-----------------------------	--	-----------


Table 2.7. Size and weight of accessories: Open Flow Restrictor OFR

- 1 – Length of the body of OFR 3 – Diameter of the inner threading
2 – Diameter of the knurled body

Description	Length of the body, mm	Diameter of the knurled body, mm	Type of the inner threading	Diameter and pitch of the inner threading, mm	Weight, g
Open Flow Restrictor OFR	19,0±0,2	11,0±0,2	M7 x 1,0	7,0/1,0	6,5±0,5

Table 2.8. Size and weight of accessories: Spray Tip Cleaning Adapter TCA

A – Spray Tip Cleaning Adapter TCA, disassembled view

B – Spray Tip Cleaning Adapter TCA, work-ready view

1 – Length of the adapter

2 – Diameter of the knurled body

3 – Diameter of the outlet

4 – Diameter of the inner threading


5 – Diameter of the body of the lock screw

6 – Length of the threading of the lock screw

7 – Diameter of the threading of the lock screw

Description	Length of the adapter, mm	Diameter of the knurled body, mm	Diameter of the outlet, mm	Type of the inner threading of the body of the adapter, mm	Diameter and pitch of the inner threading, mm	Diameter of the body of the lock screw, mm	Length/diameter/pitch of the threading of the lock screw, mm	Weight, g
Spray Tip Cleaning Adapter TCA	24,0±0,2	14,0±0,2	7,0	M12	12,0 x 1,0	14,0±0,2	6,0/14,0/1,0	15,0±1,0

Table 2.9. Size and weight of accessories: Lexan disc with four cones

1 – Diameter of the disc

2 – Height of the disc


3 – Diameter of the base of the cones

4 – Diameter of impact openings

5 – Height of the side of the disc

Description	Diameter of the disc, mm	Height of the disc, mm	Diameter of the base of the cones	Diameter of the impact openings, mm	Height of the side of the disc, mm	Weight, g
Lexan disc with four cones	115,0±0,2	25,0±0,2	35,0±0,5 37,0±0,5 42,0±0,5 44,0±0,5	10,0±0,5; 8,0±0,5; 6,0±0,5; 3,0±0,5	8,0	20,0±2,0

Table 2.10. Size and weight of accessories: Neoprene cones, a set of 6 ea

- 1 – Outer diameter of the base of a cone
 2 – Outer diameter of the blunted tip of a cone
 3 – Height of a cone

- 4 – Diameter of the opening of the base of a cone
 5 – Diameter of the impact opening

Наименование	Outer diameter of the base of a cone, mm	Outer diameter of the blunted tip of a cone, mm	Height of a cone, mm	Diameter of the opening of the base of a cone, mm	Diameter of the impact opening, mm	Weight, g
Neoprene cones, a set of 6 ea including:						
– Neoprene cone No. 1	10,0	22,0	19,0±0,5	22,0±1,0	10,0±1,0	3±0,5;
– Neoprene cone No. 2	16,0	28,0	19,0±0,5	28,0±1,0	16,0±1,0	4,5±0,5;
– Neoprene cone No. 3	22,0	37,0	25,0±0,5	37,0±1,0	22,0±1,0	9,0±0,5;
– Neoprene cone No. 4	29,0	46,0	30,0±0,5	47,0±1,0	29,0±1,0	15,0±0,5;
– Neoprene cone No. 5	38,0	58,0	34,0±0,5	58,0±1,0	38,0±1,0	25,0±0,5;
– Neoprene cone No. 6	46,0	71,0	39,0±0,5	71,0±1,0	46,0±1,0	38,0±0,5

Table 3. Parameters of threading

Name of parts	Type of threading	Screw pitch, mm	Length of threading, mm	Depth of threading, mm
The lid of the can	M50	2,0±0,05	10,0±0,05	1,5±0,05
Fittings of the closed contact probes CP SP, CP 1 – 20 mm	M7	1,0±0,05	6,0±0,05	1,0±0,05
Fittings of cervical probes: – Cervical probe 19x5 mm – Cervical probe 19x15 mm – Cervical probe 25x5 mm – Cervical probe 19 mm, dome shaped – Cervical probe 25 mm, dome shaped	M7	1,0±0,05	6,0±0,05	1,0±0,05
Fittings of open spray tips OS SS, OS BS; OS SOFT; OS LL; OS A, OS B, OS C, OS D	M7	1,0±0,05	6,0±0,05	1,0±0,05
Fittings of Open Flow Restrictor OFR and Spray Tip Cleaning Adapter TCA	M7	1,0±0,05	6,0±0,05	1,0±0,05
Fitting of the brass body of the Open Spray Tip OS LL	M7	1,0±0,05	6,0±0,05	1,0±0,05

Table 4. Parameters of knurling and roughness

Names of parts	Type of knurling	Pitch of knurling, mm	Height of knurling, mm	Roughness of the contact surface, μm
----------------	------------------	-----------------------	------------------------	--------------------------------------

Brass thimble for connecting the tips to the discharge tube of the exhaust valve	Reticular	1,5±0,1	0,5±0,1	1,5±0,1
The knurled brass bodies of: – Closed contact probes 1-20 mm – Cervical probes: – Cervical probe 19x5 mm – Cervical probe 19x15 mm – Cervical probe 25x5 mm – Cervical probe 19 mm, dome shaped – Cervical probe 25 mm, dome shaped – Open spray tip - OS SS, – Open spray tip - OS BS, – Open spray tips OS A, OS B, OS C, OS D, – Open spray tip - OS SOFT, – Open Flow Restrictor OFR	Reticular	1,5±0,1	0,5±0,1	1,5±0,1
Open Spray Tip OS LL: – Knurled brass body of the tip – Knurled stainless steel adapter body with Luer Lock	Reticular	1,5±0,1	0,5±0,1	1,5±0,1
Spray Tip Cleaning Adapter TCA: – Knurled adapter body – Knurled lock screw	Reticular	1,5±0,05	0,5±0,1	1,5±0,1
Teflon® coating of the cervical probes' heads	---	---	---	1,0±0,2

3. Indications for use of the medical device

Many types of benign skin lesions can provide indications for cryotherapy with the help of the cryosurgical instrument CryoPro: all types of warts (common warts, flat warts, plantar warts, pointed condyloma and senile warts), acne, molluscum contagiosum, seborrheic keratosis, dermatofibroma, angioma, myxoid cyst, sebaceous adenoma, pyogenic granuloma, cutaneous horn, papilloma, lentigo, and also some forms of precancerous pathologies: actinic keratosis, malignant lentigo, some types of papilloma, initial stages of basal cell carcinoma, and small single lesions suspicious for melanoma. With the help of cryosurgery, the ingrown nails and some types of hemorrhoid and cervical cancer may also be treated. If there is a suspicion about a malignant character of a skin lesion, a careful study and a consultation from an experienced dermatologist or oncologist are quite necessary. Cryotherapy of such lesions should be executed by most experienced cryosurgeons only.

4. Counter-indications

During the operational period from the beginning of the production of cryosurgical instruments CryoPro in 1986 there have been found no counter-indications as to their use. A risk of injury to

both patients and users is non-existent when the users observe directions for transportation, storage and operation of the cryosurgical instruments. There have not been found any counter-indications for the use of the cryosurgical instrument under pregnancy. The age of the patient is also not a counter-indication for the use, with the exception of children under 1 year who may react too strong to irritations. Light mechanical or thermal damages of skin are not a counter-indication for use. If more extensive damages of skin are seen, a consultation of a dermatologist is needed.

General counter-indications for cryosurgery are: individual hyper-sensitivity to cold, especially when it causes urticaria; the hyper-pigmented skin; the lesions located on the eyelid margins or on the nasolabial fold or on the wings of nose, or on the haired areas; cryoglobulinemia; collagen or autoimmune diseases; Raynaud's disease; multiple myeloma; pyoderma gangrenosum; advanced basal cell carcinoma, and any extensive damages of the skin surface.

5. Safety measures

During the operational life time of the cryosurgical instrument CryoPro in its specified manufacturing modifications one should observe the following rules of its safe operation:

- Only a doctor may execute the cryosurgical procedures in a strict compliance with the technique of the corresponding procedures and with respect to all known indications and counter-indications applicable to their use;
- Cryosurgical instrument CryoPro is intended to be used with liquid nitrogen only! Never use any different cryogenic liquids with this instrument!
- Never use the instrument contrary to its intended use specified in the User Manual;
- Store the liquid nitrogen supply solely in Dewar containers and solely at the specially assigned locations having suction-and-exhaust ventilation, in cool places off from heat sources and direct sunlight; take measures against occasional tilting of Dewar containers. The containers should be periodically cleaned from accumulating particles and ice;
- Fill the can of the cryosurgical instrument with liquid nitrogen observing all safety measures (protective gloves and glasses) and using a detachable withdrawal tube or other devices that are not part of the delivery set and supplied separately; do not splash or spill liquid nitrogen on objects that could be damaged by it. The withdrawal devices should be equipped with a working LN filter to prevent clogging of tips by foreign particles that liquid nitrogen can contain;
- Fill the can at no more than 70% of its full volume leaving free space for the liquid nitrogen vapor to create the sufficient working pressure in it (about 0,7 bar);
- After filling the can and right after the liquid nitrogen stops boiling put the lid of the can back in its place screwing it with a moderate force of the hand; never screw or unscrew

the lid using the discharge tube of the exhaust valve as a lever because it can be bent or broken;

- Never open the exhaust valve when a tip is not connected to the discharge tube; better use the Open Flow Restrictor OFR permanently (Picture 30, Page 31);
- Be enough careful working with tips in the vicinity of the eyes taking all precautions against liquid nitrogen affecting the eye tissues or accidental touching the eyes by a tip. If such procedures ever become necessary we recommend using simple medical devices that can help to protect the eyes;
- When during working hours there are children in the medical facility, the doctor should watch that they could not take the cryosurgical device and use it uncontrollably;
- During temporary breaks within a working day keep the cryosurgical instrument **vertically** (with its exhaust valve looking upwards) in a cool place off the heat sources and direct sunlight; make sure that the plastic supporting base of the can is always in its place to reduce possibility of occasional tilting of the device;
- If the device that contain liquid nitrogen fell or was overthrown, immediately get it back into the vertical position, watch its behavior for 15-20 minutes and then check its operational capability;
- If any large areas of frosting or ice appear, or are otherwise discovered, on the surface of the can and an intensive permanent hissing of the safety valve is heard you have to put on protective glasses and gloves, take the tip off the discharge tube and, pressing the trigger of the exhaust valve, let the excess pressure come out the can, then open the can lid and remove any remaining liquid nitrogen from it. Such device has to be necessarily sent for servicing to the manufacturer through its authorized representative (the address is shown on Page 43). Light hissing of the safety valve without frosting or ice on the can's surface does not report any malfunction;
- In case any outer mechanical damages of the device's components appeared or are otherwise discovered (for example, deformed or broken parts of the exhaust valve, safety valve, the lid of the can, or deformed/pierced outer body of the can), stop using the cryosurgical instrument and send it for servicing or repair to the manufacturer through its authorized representative (the address is shown on Page 43);
- In case of clogging of spray tips by foreign particles they have to be cleaned-through using the Spray Tip Cleaning Adapter TCA (Picture 34, Page 32). Do not clean them using occasional and unreliable stuff (for example, wooden or plastic sticks or copper wire);
- If the room in which people work with the cryosurgical instrument does not have sufficient ventilation (for example, in winter) it is recommended to air it periodically during the working hours to make sure there is no excessive concentration of nitrogen in the environment.
- At the end of the working hours decant any remaining liquid nitrogen from the device can into a Dewar container;

the lid using the discharge tube of the exhaust valve as a lever because it can be bent or broken;

- Never open the exhaust valve when a tip is not connected to the discharge tube; better use the Open Flow Restrictor OFR permanently (Picture 30, Page 31);
- Be enough careful working with tips in the vicinity of the eyes taking all precautions against liquid nitrogen affecting the eye tissues or accidental touching the eyes by a tip. If such procedures ever become necessary we recommend using simple medical devices that can help to protect the eyes;
- When during working hours there are children in the medical facility, the doctor should watch that they could not take the cryosurgical device and use it uncontrollably;
- During temporary breaks within a working day keep the cryosurgical instrument **vertically** (with its exhaust valve looking upwards) in a cool place off the heat sources and direct sunlight; make sure that the plastic supporting base of the can is always in its place to reduce possibility of occasional tilting of the device;
- If the device that contain liquid nitrogen fell or was overthrown, immediately get it back into the vertical position, watch its behavior for 15-20 minutes and then check its operational capability;
- If any large areas of frosting or ice appear, or are otherwise discovered, on the surface of the can and an intensive permanent hissing of the safety valve is heard you have to put on protective glasses and gloves, take the tip off the discharge tube and, pressing the trigger of the exhaust valve, let the excess pressure come out the can, then open the can lid and remove any remaining liquid nitrogen from it. Such device has to be necessarily sent for servicing to the manufacturer through its authorized representative (the address is shown on Page 43). Light hissing of the safety valve without frosting or ice on the can's surface does not report any malfunction;
- In case any outer mechanical damages of the device's components appeared or are otherwise discovered (for example, deformed or broken parts of the exhaust valve, safety valve, the lid of the can, or deformed/pierced outer body of the can), stop using the cryosurgical instrument and send it for servicing or repair to the manufacturer through its authorized representative (the address is shown on Page 43);
- In case of clogging of spray tips by foreign particles they have to be cleaned-through using the Spray Tip Cleaning Adapter TCA (Picture 34, Page 32). Do not clean them using occasional and unreliable stuff (for example, wooden or plastic sticks or copper wire);
- If the room in which people work with the cryosurgical instrument does not have sufficient ventilation (for example, in winter) it is recommended to air it periodically during the working hours to make sure there is no excessive concentration of nitrogen in the environment.
- At the end of the working hours decant any remaining liquid nitrogen from the device can into a Dewar container;

- If more than one person is working with the same cryosurgical instrument they have to inform each other about any noticed malfunctions/defects of the device.

6. Pre-starting and operating procedures

6.1. General directions

Unscrew the lid of the can by turning it anti-clockwise. Place the open can of the CryoPro under the filter nozzle of the withdrawal tube and fill the can up to no more than 70% of its volume. Operational procedures for the LN withdrawal devices may be different depending on their model. Prior to using them, you should read the corresponding manual instructions.

After filling the can of the cryosurgical instrument CryoPro with liquid nitrogen and installing the lid back in its place you can start doing cryosurgical procedures. To proceed, a certain tip should be connected to the discharge tube of the cryosurgical instrument. Doctor chooses both type and size of the tips depending on the character of a skin lesion to be treated. All contact and cervical probes (Picture 6-21, Page 21-28) and open spray tips (Picture 22-29, Page 29-30) have threaded fittings and connect to the cryosurgical instrument CryoPro with the help of the knurled thimble by the force of the fingers. There should be no excessive force applied while doing that. The period of static storage of liquid nitrogen in the can of the device (until the complete evaporation with its lid closed tightly and with no consumption for procedures) is: for CryoPro Maxi – 20-24 hours, for CryoPro Mini – 12-14 hours.

6.2. Lesion freezing procedure


The freezing of skin lesions is performed using closed contact probes and open spray tips. Contact and cervical probes (Picture 6-21, Page 21-28) affect the procedure area by a direct contact between the cooled probe's head and skin, and give a deeper freezing with its minimal lateral spread. The size of a probe is to be chosen depending on the size of a lesion. The lesion should be moistened with a drop of water or gel before the procedure to increase the heat exchange between the skin surface and the probe's head. Place the probe's head onto the skin when it is still warm, press it slightly and push the trigger of the exhaust valve. When going through the probe the liquid nitrogen is cooling the head from inside and then comes out into the air through a special silicone tube (waste tube).

The silicone tube that is attached to a probe is intended to discharge the cold gas in a direction away from the patient, and it should be given a safe direction from the very beginning of a procedure. The tube can be easily cut when necessary. In case the trigger of the exhaust valve is stuck in an open position (it is very unlikely though) you should immediately let the excess pressure out from the can slowly unscrewing its lid. This procedure should be done while the head of a probe is frozen to the treatment area.

Most lesions could be sufficiently frozen in less than 45 seconds. Yet freezing of the mucous membranes should never last longer than 60 seconds.

Give approximately 5 seconds for the probe to get well cooled and keep freezing for about one minute while an ice ring of about 1 mm wide is formed around. Then stop the freezing and wait few seconds to ensure the head could be easily separated from skin. The procedure using contact probes and "fast freezing – slow thawing" cycles ensures a very efficient treatment with a high rate of impact on pathological tissues and great accuracy.


6.3. Closed contact probes (Picture 6-16, Pages 25-30)

Picture 6. Closed contact probe - CP SP

Intended use

The closed contact probe - CP SP is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on small lesion with a cavity in the center – for example, on angiomas or on molluscum contagiosum. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating

Picture 7. Closed contact probe - CP 1 mm.


Intended use

The closed contact probe - CP 1 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on very small warts, angiomas, lesions around the eyes, and trichiasis. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an

Most lesions could be sufficiently frozen in less than 45 seconds. Yet freezing of the mucous membranes should never last longer than 60 seconds.

Give approximately 5 seconds for the probe to get well cooled and keep freezing for about one minute while an ice ring of about 1 mm wide is formed around. Then stop the freezing and wait few seconds to ensure the head could be easily separated from skin. The procedure using contact probes and "fast freezing – slow thawing" cycles ensures a very efficient treatment with a high rate of impact on pathological tissues and great accuracy.


6.3. Closed contact probes (Picture 6-16, Pages 25-30)

Picture 6. Closed contact probe - CP SP

Intended use

The closed contact probe - CP SP is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on small lesion with a cavity in the center – for example, on angiomas or on molluscum contagiosum. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating

Picture 7. Closed contact probe - CP 1 mm.


Intended use

The closed contact probe - CP 1 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on very small warts, angiomas, lesions around the eyes, and trichiasis. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an

Most lesions could be sufficiently frozen in less than 45 seconds. Yet freezing of the mucous membranes should never last longer than 60 seconds.

Give approximately 5 seconds for the probe to get well cooled and keep freezing for about one minute while an ice ring of about 1 mm wide is formed around. Then stop the freezing and wait few seconds to ensure the head could be easily separated from skin. The procedure using contact probes and "fast freezing – slow thawing" cycles ensures a very efficient treatment with a high rate of impact on pathological tissues and great accuracy.


6.3. Closed contact probes (Picture 6-16, Pages 25-30)

Picture 6. Closed contact probe - CP SP

Intended use

The closed contact probe - CP SP is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on small lesion with a cavity in the center – for example, on angiomas or on molluscum contagiosum. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating


Picture 7. Closed contact probe - CP 1 mm.

Intended use

The closed contact probe - CP 1 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on very small warts, angiomas, lesions around the eyes, and trichiasis. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an

with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.


Picture 8. Closed contact probe – CP 2 mm.

Intended use

The closed contact probe – CP 2 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on small flat warts, lesions around the eyes, and trichiasis. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.


inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 9. Closed contact probe - CP 3 mm.

Intended use


The closed contact probe - CP 3 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 10. Closed contact probe – CP 4 mm.

Intended use


The closed contact probe – CP 4 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 11. Closed contact probe – CP 5 mm.

Intended use


The closed contact probe – CP 5 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 10. Closed contact probe – CP 4 mm.

Intended use


The closed contact probe – CP 4 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 11. Closed contact probe – CP 5 mm.

Intended use


The closed contact probe – CP 5 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 10. Closed contact probe – CP 4 mm.

Intended use


The closed contact probe – CP 4 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 11. Closed contact probe – CP 5 mm.

Intended use


The closed contact probe – CP 5 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 10. Closed contact probe – CP 4 mm.

Intended use


The closed contact probe – CP 4 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 11. Closed contact probe – CP 5 mm.

Intended use


The closed contact probe – CP 5 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 12. Closed contact probe – CP 6 mm.

Intended use


The closed contact probe – CP 6 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 13. Closed contact probe – CP 8 mm.

Intended use


The closed contact probe – CP 8 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on flat warts, lesions around the eyes, and lentiginosis if their size is similar to the size of the probe's head. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 14. Closed contact probe – CP 10 mm.

Intended use

The closed contact probe – CP 10 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use with some types of basal cell carcinoma and hemorrhoid. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

Picture 15. Closed contact probe – CP 15 mm.

Intended use

The closed contact probe – CP 15 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on sizeable basal cell carcinomas and hemorrhoid. It may also be used for different indications that are matching its properties and configuration, and the practical use of this type of a probe. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.


Picture 16. Closed contact probe – CP 20 mm.

Intended use


The closed contact probe – CP 20 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on sizeable basal cell carcinomas. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

6.4. Cervical probes


Cervical probes of five sizes (Pictures 17-21, Page 26-28) are used for treatment of lesions located on mucous membranes. Their use is similar to other contact probes. The timing of a procedure is set by an experienced doctor depending on the character of a lesion.

Picture 17. Cervical probe 19x5 mm.

Picture 18. Cervical probe 19x15 mm.


Picture 16. Closed contact probe – CP 20 mm.

Intended use


The closed contact probe – CP 20 mm is intended to be used for freezing of skin lesions while under their direct contact with the tip's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for use on sizeable basal cell carcinomas. The tip's head has a smooth galvanic coating with gold that is an inert non-toxic material; it gets detached from skin pretty easily during the thaw phase.

6.4. Cervical probes

Cervical probes of five sizes (Pictures 17-21, Page 26-28) are used for treatment of lesions located on mucous membranes. Their use is similar to other contact probes. The timing of a procedure is set by an experienced doctor depending on the character of a lesion.

Picture 17. Cervical probe 19x5 mm.

Picture 18. Cervical probe 19x15 mm.

Intended use

The cervical probe 19x5 mm is intended to be used for freezing of tissue while under its direct contact with the probe's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for removal of lesions located on mucous membranes by using the ultra-low temperature of liquid nitrogen if there are valid indications for that type of treatment and if the diameter and the depth of pathology match the probe's configuration. The tip's head has a smooth protective Teflon coating that prevents damage to the healthy skin and is easily detachable from the tender cervical mucosa during the thaw phase.

Picture 19. Cervical probe 25x5 mm.

Intended use

The cervical probe 25x5 mm is intended to be used for freezing of tissue while under its direct contact with the probe's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for removal of lesions located on mucous membranes by using the ultra-low temperature of liquid nitrogen if there are valid indications for that type of treatment and if the diameter and the depth of pathology match the probe's configuration. The tip's head has a smooth protective Teflon coating that prevents damage to the healthy skin and is easily detachable from the tender cervical mucosa during the thaw phase.

Intended use

The cervical probe 19x15 mm is intended to be used for freezing of tissue while under its direct contact with the probe's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for removal of lesions located on mucous membranes by using the ultra-low temperature of liquid nitrogen if there are valid indications for that type of treatment and if the diameter and the depth of pathology match the probe's configuration. The tip's head has a smooth protective Teflon coating that prevents damage to the healthy skin and is easily detachable from the tender cervical mucosa during the thaw phase.

Picture 20. Cervical probe 19 mm, dome shaped.

Intended use

The cervical probe 19 mm, dome shaped is intended to be used for freezing of tissue while under its direct contact with the probe's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for removal of lesions located on mucous membranes by using the ultra-low temperature of liquid nitrogen if there are valid indications for that type of treatment and if the diameter and the depth of pathology match the probe's configuration. The tip's head has a smooth protective Teflon coating that prevents damage to the healthy skin and is easily detachable from the tender cervical mucosa during the thaw phase.

Picture 21. Cervical probe 25 mm, dome shaped.


Intended use

The cervical probe 25 mm, dome shaped is intended to be used for freezing of tissue while under its direct contact with the probe's head that is cooled by liquid nitrogen circulating inside the probe. It is recommended for removal of lesions located on mucous membranes by using the ultra-low temperature of liquid nitrogen if there are valid indications for that type of treatment and if the diameter and the depth of pathology match the probe's configuration. The tip's head has a smooth protective Teflon coating that prevents damage to the healthy skin and is easily detachable from the tender cervical mucosa during the thaw phase.

6.5. Open spray tips

Freezing of lesions with the help of open spray tips (Picture 22-29, Page 29-30) is done by a direct impact of a subtle spray of liquid nitrogen streaming out the orifice of the tip. Doctor directs the tip at the procedure area and presses the trigger of the exhaust valve. Freezing time using the open spray can be from 5 to 20 seconds. For a procedure to produce a better effect, it is recommended to perform several "freeze-thaw" cycles. Open spray tips - OS A and OS B which have the diameter of the orifice of 0,8 and 0,6 mm respectively allow treating relatively large and thick lesions and could also be used for desquamation. Open spray tip - OS C with the diameter of its orifice of 0,45 mm is recommended for treatment of small warts and keratolytic lesions. Open spray tip - OS D with the diameter of its orifice of 0,3 mm is used for treatment of smallest lesions including warts and small papilloma. Open spray tip - OS SS (straight) is used where it is by some reason impossible to use open tips OS A, B, C or D – for example, inside the mouth. Open spray tip - OS BS has a bent ending and allows to treat lesions in relatively inconvenient places – on the nose or in the helix etc.


Open spray tip - OS SOFT is used for treating acne and for soft cryo-peeling. Open spray tip OS LL is used for the syringe needles to be connected to the cryosurgical instrument if it is ever needed to inject liquid nitrogen into voluminous lesions.

Picture 22. Open spray tip – OS SS.

Intended use


The open spray tip – OS SS is intended to treat lesions located in natural hollows or orifices of the human body, i.e. in such places where it is difficult or impossible to use aperture tips. A thin open spray of liquid nitrogen is applied direct to skin without any contact between the tip and the skin. The tip is 86 mm long and its orifice is 0,4 mm.

Picture 23. Open spray tip - OS BS.

Intended use


The open spray tip – OS BS (bent) is intended to treat lesions located in natural hollows or orifices of the human body that are inaccessible for the straight open spray tip. A thin open spray of liquid nitrogen is applied direct to skin without any contact between the tip and the skin. The tip is 83 mm long and its orifice is 0,4 mm.

Picture 24. Open spray tip - OS SOFT.

Intended use


The open spray tip – OS SOFT is intended to help patients with akne. It has a row of nine orifices with the diameter of 0,22 mm each and produces an open flat spray of liquid nitrogen of high dispersity and soft on impact. Such a spray allows to make a softest peeling of the very surface of skin. A thin open spray of liquid nitrogen is applied direct to skin without any contact between the tip and the skin.

Picture 25. Open spray tip - OS A.

Intended use


The open spray tip - OS A (aperture) is intended to treat skin lesions by the open spray of liquid nitrogen without contact with patient's skin. With it, one can treat any skin lesions that are located on easily accessible parts of human body. Aperture OS A has inner diameter of 0,8 mm.

Picture 26. Open spray tip - OS B.

Intended use


The open spray tip - OS B (aperture) is intended to treat skin lesions by the open spray of liquid nitrogen without contact with patient's skin. With it, one can treat any skin lesions that are located on easily accessible parts of human body. Aperture OS B has inner diameter of 0,6 mm.


Picture 27. Open spray tip - OS C.

Intended use

The open spray tip - OS C (aperture) is intended to treat skin lesions by the open spray of liquid nitrogen without contact with patient's skin. With it, one can treat any skin lesions that are located on easily accessible parts of human body. Aperture OS C has inner diameter of 0,45 mm.


Picture 28. Open spray tip - OS D.

Picture 29. Open spray tip - OS LL.

Intended use


The open spray tip - OS LL allows connecting the syringe needles to cryosurgical instruments CryoPro Maxi and CryoPro Mini. It consists of two parts: brass tip and stainless steel Luer Lock adapter. Brass tip connects to the discharge tube of the CryoPro and needles connect to the Luer Lock adapter. That can be necessary in cases when the liquid nitrogen spray should be injected into the lesions. The tip does not contact patient's skin.

Picture 26. Open spray tip - OS B.

Intended use


The open spray tip - OS B (aperture) is intended to treat skin lesions by the open spray of liquid nitrogen without contact with patient's skin. With it, one can treat any skin lesions that are located on easily accessible parts of human body. Aperture OS B has inner diameter of 0,6 mm.

Picture 27. Open spray tip - OS C.

Intended use


The open spray tip - OS C (aperture) is intended to treat skin lesions by the open spray of liquid nitrogen without contact with patient's skin. With it, one can treat any skin lesions that are located on easily accessible parts of human body. Aperture OS C has inner diameter of 0,45 mm.

Picture 28. Open spray tip - OS D.

Intended use

The open spray tip - OS D (aperture) is intended to treat skin lesions by the open spray of liquid nitrogen without contact with patient's skin. With it, one can treat any skin lesions that are located on easily accessible parts of human body. Aperture OS D has inner diameter of 0,3 mm.

Picture 29. Open spray tip - OS LL.

Intended use


The open spray tip - OS LL allows connecting the syringe needles to cryosurgical instruments CryoPro Maxi and CryoPro Mini. It consists of two parts: brass tip and stainless steel Luer Lock adapter. Brass tip connects to the discharge tube of the CryoPro and needles connect to the Luer Lock adapter. That can be necessary in cases when the liquid nitrogen spray should be injected into the lesions. The tip does not contact patient's skin.

The most effective freezing with open spray tips is done when they are held pretty close to the skin surface (at a distance of 5 to 10 mm) because the temperature of spray is always lower at a closer distance. It is recommended to treat the prominent lesions at a tangent to them and from two directions minimum. This technique allows to form an ice ball deeper into the lesion with a minimal impact on the adjacent tissue. A treatment with spray at 90° to skin causes a greater lateral than deeper freezing. To avoid that, apply spray on small areas interruptedly and use a tip of smaller size.

Desquamation of surface can be done using tips with larger orifice like OS A and OS B or OS SOFT from a distance of 5 cm (at that distance the spray contains more vapor) "painting" the surface with evaporating liquid nitrogen.

6.6. Additional tips (adapters)


Open Flow Restrictor OFR (Picture 30, Page 31) is mounted onto the discharge tube of the cryosurgical instrument (Picture 31-33, Page 32) and is designed to provide more safety. It is preventing incidental discharges of liquid nitrogen when there is no tip mounted on the device. Cryosurgical device cannot be used if any tip is not connected to it. It could be useful for preventing a strong discharge of liquid nitrogen at any incidental trigger of the exhaust valve when there is no tip on the discharge tube.


Picture 30. Open Flow Restrictor OFR.

Intended use


The Open Flow Restrictor OFR prevents using the cryosurgical instrument without a tip or a probe. It can be permanently connected to the discharge tube of the cryosurgical instrument. When doctor presses the trigger the liquid nitrogen will come out only when a tip or probe is connected to the OFR.

Picture 31.


Picture 32.


Picture 33.


Spray Tip Cleaning Adapter TCA (Picture 34, Page 32) is used for cleaning of open spray tips OS A, B, C and D in case they get clogged. It has a design creating an increased pressure within its body. A clogged tip is inserted into its lock screw (Picture 37, Page 33), then its body is connected to the discharge tube of the cryosurgical instrument, its lock screw is screwed on (Picture 38, Page 33) and the blowing through is being done (Picture 42, Page 34).


Picture 34. Spray Tip Cleaning Adapter TCA.

Picture 31.


Picture 32.

Picture 33.

Spray Tip Cleaning Adapter TCA (Picture 34, Page 32) is used for cleaning of open spray tips OS A, B, C and D in case they get clogged. It has a design creating an increased pressure within its body. A clogged tip is inserted into its lock screw (Picture 37, Page 33), then its body is connected to the discharge tube of the cryosurgical instrument, its lock screw is screwed on (Picture 38, Page 33) and the blowing through is being done (Picture 42, Page 34).

Picture 34. Spray Tip Cleaning Adapter TCA.

A – disassembled adapter,
B – ready-to-work adapter.

Intended use

The Spray Tip Cleaning Adapter TCA is intended to blow-through the open spray tips type OS A, B, C and D in case of their clogging. The spray tips are inserted into the lock screw, then it is screwed into the adapter body and then the adapter is connected to the exhaust pipe of the cryosurgical instrument, the trigger of the exhaust valve is pressed and the blowing by the flow of liquid nitrogen proceeds.


Picture 35.


Picture 36.

Picture 37.


Picture 38.


Picture 40.

Picture 39.

Picture 41.


Picture 42.

6.7. Protective accessories

When doing freezing procedures, especially close to eyes or where any damage to healthy skin is inappropriate, it is recommended to use protective accessories – Lexan disc and Neoprene cones. **Lexan disc with four cones** (Picture 43, Page 34) is intended for protecting the surface of healthy skin when the spray of liquid nitrogen is applied to skin. The disc is made of durable heat resistant plastic Lexan™. It has four hollows in the form of blunted cones having openings of different diameter. During a procedure the chosen cone is pressed to skin and its opening embraces the lesion which is then treated with the liquid nitrogen spray. In this case the healthy skin is protected from any damage by ultra-low temperatures..


Picture 43. Lexan disc with four cones.

Picture 44. Using lexan disc.

Neoprene cones (Picture 45, Page 35) (a set of six pieces of different sizes) realize the same function. They are made of neoprene rubber and are flexible, and therefore their form can be adjusted to match the treatment area. Lexan disc and Neoprene cones are used in procedures that employ open spray tips to remove any damage to healthy skin surrounding the treatment area.


Picture 45. Neoprene cones, a set of 6.

- Neoprene cone No. 1, large opening 22 mm, impact opening 10 mm – 1 ea
- Neoprene cone No. 2, large opening 28 mm, impact opening 16 mm – 1 ea
- Neoprene cone No. 3, large opening 37 mm, impact opening 22 mm – 1 ea
- Neoprene cone No. 4, large opening 47 mm, impact opening 29 mm – 1 ea
- Neoprene cone No. 5, large opening 58 mm, impact opening 38 mm – 1 ea
- Neoprene cone No. 6, large opening 71 mm, impact opening 46 mm – 1 ea.


Picture 44. Using lexan disc.

Neoprene cones (Picture 45, Page 35) (a set of six pieces of different sizes) realize the same function. They are made of neoprene rubber and are flexible, and therefore their form can be adjusted to match the treatment area. Lexan disc and Neoprene cones are used in procedures that employ open spray tips to remove any damage to healthy skin surrounding the treatment area.

Picture 45. Neoprene cones, a set of 6.

- Neoprene cone No. 1, large opening 22 mm, impact opening 10 mm – 1 ea
- Neoprene cone No. 2, large opening 28 mm, impact opening 16 mm – 1 ea
- Neoprene cone No. 3, large opening 37 mm, impact opening 22 mm – 1 ea
- Neoprene cone No. 4, large opening 47 mm, impact opening 29 mm – 1 ea
- Neoprene cone No. 5, large opening 58 mm, impact opening 38 mm – 1 ea
- Neoprene cone No. 6, large opening 71 mm, impact opening 46 mm – 1 ea.

Picture 46. Using a cone.

7. Storage precautions for the cryosurgical instrument

It is strongly recommended that at the end of the working day the cryosurgical instrument CryoPro be kept with its lid closed to prevent moisture condensation inside the can and ingress of any foreign particles in it. Accumulation of condensate may lead to blocking of valve mechanisms.

Cryosurgical instrument CryoPro has to be kept in dry rooms and should not experience long exposure to dust, aggressive gases, and acid or alkaline vapors. The following environmental conditions should be observed during the storage:

- Ambient temperature from -10°C to $+60^{\circ}\text{C}$
- Relative humidity from 30% to 95% under $+25^{\circ}\text{C}$
- Atmospheric pressure from 84 to 106,7 kPa.

8. Maintenance, disinfection and disposal of the cryosurgical instrument

Open spray tips do not require sterilization because they do not enter into contact with patient skin. Nevertheless, if their sterilization is necessary you should use the same procedure that is used for the closed contact probes and cervical probes.

Heads of the closed contact probes and cervical probes are contacting the patient skin directly, therefore they must be always sterilized after every use. Recommended method of sterilization is the steam autoclaving during 15 minutes at $+121^{\circ}\text{C}$.

To clean the body and outer parts of the cryosurgical instrument, the disinfecting kleenex or sprays could be first used, than a clean cloth to wipe them dry.

When the operating period is ending or in case of complete incapacity for use, the cryosurgical instrument CryoPro should be disposed of according to the rules of utilization valid for such devices on the territory of Russian Federation.

9. Transportation of the cryosurgical instrument

The cryosurgical instrument in its specified manufacturing modifications could be transported by all means of transport, in closed vehicles in accordance with valid rules of transportation, under temperatures from -20°C to +50°C and relative humidity from 40% to 95% under +25°C.

10. Complete delivery set of the cryosurgical instrument

10.1. Cryosurgical instrument CryoPro, manufacturing modification CryoPro Maxi:

1. Cryosurgical instrument CryoPro Maxi – 1 ea
2. A set of open spray tips including:
 - Open spray tip - OS BS – 1 ea
 - Open spray tip - OS A – 1 ea
 - Open spray tip - OS B – 1 ea
 - Open spray tip - OS C – 2 ea
 - Open spray tip - OS D – 1 ea
 - Stand for spray tips – 1 ea.
3. Accessories for the cryosurgical instrument CryoPro, manufacturing modification CryoPro Maxi:
 - 1) Closed contact probe - CP SP – 1 ea
 - 2) Closed contact probe - CP 1 mm – 1 ea
 - 3) Closed contact probe - CP 2 mm – 1 ea
 - 4) Closed contact probe - CP 3 mm – 1 ea
 - 5) Closed contact probe - CP 4 mm – 1 ea
 - 6) Closed contact probe - CP 5 mm – 1 ea
 - 7) Closed contact probe - CP 6 mm – 1 ea
 - 8) Closed contact probe - CP 8 mm – 1 ea
 - 9) Closed contact probe - CP 10 mm – 1 ea
 - 10) Closed contact probe - CP 15 mm – 1 ea
 - 11) Closed contact probe - CP 20 mm – 1 ea
 - 12) Cervical probe 19x5 mm – 1 ea
 - 13) Cervical probe 19x15 mm – 1 ea
 - 14) Cervical probe 25x5 mm – 1 ea
 - 15) Cervical probe 19 mm, dome shaped – 1 ea

16) Cervical probe 25 mm, dome shaped	- 1 ea
17) Open spray tip - OS SS	- 1 ea
18) Open spray tip - OS BS	- 1 ea
19) Open spray tip - OS SOFT	- 1 ea
20) Open spray tip - OS A	- 1 ea
21) Open spray tip - OS B	- 1 ea
22) Open spray tip - OS C	- 1 ea
23) Open spray tip - OS D	- 1 ea
24) Open spray tip - OS LL	- 1 ea
25) Open Flow Restrictor OFR	- 1 ea
26) Spray Tip Cleaning Adapter TCA	- 1 ea
27) Lexan disc with 4 cones	- 1 ea
28) Neoprene cones, a set of 6 including:	- 1 set
- neoprene cone No. 1, diameter 10 mm	- 1 ea
- neoprene cone No. 2, diameter 16 mm	- 1 ea
- neoprene cone No. 3, diameter 22 mm	- 1 ea
- neoprene cone No. 4, diameter 29 mm	- 1 ea
- neoprene cone No. 5, diameter 38 mm	- 1 ea
- neoprene cone No. 6, diameter 46 mm	- 1 ea
29) User Manual	- 1 ea

10.2. Cryosurgical instrument CryoPro, manufacturing modification CryoPro Mini:

1. Cryosurgical instrument CryoPro Mini – 1 ea
2. A set of open spray tips including:

– Open spray tip - OS BS	– 1 ea
– Open spray tip - OS A	– 1 ea
– Open spray tip - OS B	– 1 ea
– Open spray tip - OS C	– 2 ea
– Open spray tip - OS D	– 1 ea
– Stand for spray tips	– 1 ea.
3. Accessories for the cryosurgical instrument CryoPro, manufacturing modification CryoPro Mini:

1) Closed contact probe - CP SP	- 1 ea
2) Closed contact probe - CP 1 mm	- 1 ea
3) Closed contact probe - CP 2 mm	- 1 ea
4) Closed contact probe - CP 3 mm	- 1 ea
5) Closed contact probe - CP 4 mm	- 1 ea
6) Closed contact probe - CP 5 mm	- 1 ea
7) Closed contact probe - CP 6 mm	- 1 ea

- | | |
|---|---------|
| 8) Closed contact probe - CP 8 mm | - 1 ea |
| 9) Closed contact probe - CP 10 mm | - 1 ea |
| 10) Closed contact probe - CP 15 mm | - 1 ea |
| 11) Closed contact probe - CP 20 mm | - 1 ea |
| 12) Cervical probe 19x5 mm | - 1 ea |
| 13) Cervical probe 19x15 mm | - 1 ea |
| 14) Cervical probe 25x5 mm | - 1 ea |
| 15) Cervical probe 19 mm, dome shaped | - 1 ea |
| 16) Cervical probe 25 mm, dome shaped | - 1 ea |
| 17) Open spray tip - OS SS | - 1 ea |
| 18) Open spray tip - OS BS | - 1 ea |
| 19) Open spray tip - OS SOFT | - 1 ea |
| 20) Open spray tip - OS A | - 1 ea |
| 21) Open spray tip - OS B | - 1 ea |
| 22) Open spray tip - OS C | - 1 ea |
| 23) Open spray tip - OS D | - 1 ea |
| 24) Open spray tip - OS LL | - 1 ea |
| 25) Open Flow Restrictor OFR | - 1 ea |
| 26) Spray Tip Cleaning Adapter TCA | - 1 ea |
| 27) Lexan disc with 4 cones | - 1 ea |
| 28) Neoprene cones, a set of 6 including: | - 1 set |
| – neoprene cone No. 1, diameter 10 mm | – 1 ea |
| – neoprene cone No. 2, diameter 16 mm | – 1 ea |
| – neoprene cone No. 3, diameter 22 mm | – 1 ea |
| – neoprene cone No. 4, diameter 29 mm | – 1 ea |
| – neoprene cone No. 5, diameter 38 mm | – 1 ea |
| – neoprene cone No. 6, diameter 46 mm | – 1 ea |
| 29) User Manual | - 1 ea. |

11. Packing and marking of the cryosurgical instrument

Cryosurgical instrument in different manufacturing modifications is delivered in original packing boxes that must not have any traces of damage or opening. If any bents, breaks or other defects are discovered the equipment must not be used.

The following information is indicated on the marking label of the device (Picture 48, Page 41):

- Name of the device
- Serial number (it is on the bottom of the device, see Picture 47, Page 40)
- Text warning about the requirement to be used with liquid nitrogen only
- CE certification mark

8) Closed contact probe - CP 8 mm	- 1 ea
9) Closed contact probe - CP 10 mm	- 1 ea
10) Closed contact probe - CP 15 mm	- 1 ea
11) Closed contact probe - CP 20 mm	- 1 ea
12) Cervical probe 19x5 mm	- 1 ea
13) Cervical probe 19x15 mm	- 1 ea
14) Cervical probe 25x5 mm	- 1 ea
15) Cervical probe 19 mm, dome shaped	- 1 ea
16) Cervical probe 25 mm, dome shaped	- 1 ea
17) Open spray tip - OS SS	- 1 ea
18) Open spray tip - OS BS	- 1 ea
19) Open spray tip - OS SOFT	- 1 ea
20) Open spray tip - OS A	- 1 ea
21) Open spray tip - OS B	- 1 ea
22) Open spray tip - OS C	- 1 ea
23) Open spray tip - OS D	- 1 ea
24) Open spray tip - OS LL	- 1 ea
25) Open Flow Restrictor OFR	- 1 ea
26) Spray Tip Cleaning Adapter TCA	- 1 ea
27) Lexan disc with 4 cones	- 1 ea
28) Neoprene cones, a set of 6 including:	- 1 set
– neoprene cone No. 1, diameter 10 mm	– 1 ea
– neoprene cone No. 2, diameter 16 mm	– 1 ea
– neoprene cone No. 3, diameter 22 mm	– 1 ea
– neoprene cone No. 4, diameter 29 mm	– 1 ea
– neoprene cone No. 5, diameter 38 mm	– 1 ea
– neoprene cone No. 6, diameter 46 mm	– 1 ea
29) User Manual	- 1 ea.

11. Packing and marking of the cryosurgical instrument

Cryosurgical instrument in different manufacturing modifications is delivered in original packing boxes that must not have any traces of damage or opening. If any bents, breaks or other defects are discovered the equipment must not be used.

The following information is indicated on the marking label of the device (Picture 48, Page 41):

- Name of the device
- Serial number (it is on the bottom of the device, see Picture 47, Page 40)
- Text warning about the requirement to be used with liquid nitrogen only
- CE certification mark


- Sterilization mark (does not used by the manufacturer)
- Manufacture date (is specified by the manufacturer on the label with serial number located on the bottom of the device, see Picture 47, Page 40)
- Life time (is not indicated on the device. Specified in the User Manual)
- Manufacturing company, address of manufacturing facility, contact information.

The following is specified on the label with serial number (Picture 47, Page 40):

- Cryosurgical instrument's denomination (CryoPro Maxi or CryoPro Mini)
- Manufacturing modification of the cryosurgical instrument: L – CryoPro Maxi, S – CryoPro Mini
- Year of manufacture: the first two figures of the number
- Sequence number of the cryosurgical instrument: the next four figures of the number – they are the manufacturing sequence number of a cryosurgical instrument from the start of its production.

Picture 47. Serial numbers of cryosurgical instruments CryoPro Maxi (left) and CryoPro Mini (right)

Picture 48. Manufacturer marking labels on the cryosurgical instruments CryoPro Maxi (left) and CryoPro Mini (right).

The information specified on the marking labels is listed on Page 44.

Marking on the packing box


There is the following information on the packing boxes:

- Name of the cryosurgical instrument
- Name and address of the manufacturer
- Check-marking of the manufacturing modification (CryoPro Maxi or CryoPro Mini)
- The symbol "Read the User Manual".

Before packing, the device is cleaned, dried and covered with a thin layer of oil. The packing is made using marked cardboard boxes (Picture 49-50, Page 42).

Picture 49.

Picture 50. Packing cardboard box of the cryosurgical instrument.

12. Operational lifetime of the cryosurgical instrument

Operational Lifetime of the cryosurgical instrument is 5 years from the date of manufacture.

13. Warranty

13.1. Conditions of use of warranty

The manufacturer guarantees the quality of the cryosurgical instrument CryoPro for a period of 12 months from putting it into operation but not more than 18 months from ship date if it is used in accordance with the User Manual. The guarantee covers defects of parts and assembling. The guarantee does not cover any defects originated from failure to follow the rules of safe handling by users. Servicing is done based on purchase documents.

The guarantee does not cover defects, malfunctions and damages happened out of the following:

- Failure by a user to comply with, or negligent observation of, the rules of working with the cryosurgical instrument written in the User Manual that brought about a breakdown or failure of components of the device;
- Any user intervention into the design of the cryosurgical instrument;
- Any independent servicing or repair done by a user;
- Damages have been caused due to improper storage, packaging or transportation by user.

14. Requirements on technical servicing and repair of the medical device CryoPro

Cryosurgical instrument CryoPro does not contain any parts or components that may be repaired or replaced by users. In case the cryosurgical instrument develops a malfunction preventing its safe operation, it is necessary to contact the technical team of the authorized representative of Cortex Technology ApS in Russian Federation, MAXmedical, LLC at its address:

6 Sirenevaia Street, Office 54
Moskovskaia oblast, Schyolkovo
141109 Russian Federation.
Phone: +7 499 713 35 55.
E-mail: support@maxmedical.ru

15. Manufacturer of the medical device CryoPro

Cortex Technology ApS, Smedevaenget 10, 9560 Hadsund Denmark. Phone +45 9857 4100, fax +45 9857 2223. E-mail: cortex@cortex.dk .


16. Authorized representative of the manufacturer

Limited Liability Company "MAXmedical", 6 Sirenevaia Street, Office 54, Moskovskaia oblast, Shchyolkovo, 141109 Russian Federation. Phone +7 499 502 79 44, fax +7 499 502 50 36.
Email: info@maxmedical.ru .

If any problems with the cryosurgical instrument CryoPro occur within the warranty period please tell the model name, its serial number, the date of purchase and the description of the failure to the authorized representative of Cortex Technology ApS in Russia, MAXmedical, LLC company at the following address:

6 Sirenevaia Street, Office 54, Moskovskaia oblast, Shchyolkovo, 141109 Russian Federation.
Phone +7 499 502 79 44, +7 499 713 35 55, fax +7 499 502 50 36. Email: info@maxmedical.ru .

Cortex Technology ApS, Denmark

Gunnar B. Svendsen, CEO and Owner

25 May 2015


16. Authorized representative of the manufacturer

Limited Liability Company "MAXmedical", 6 Sirenevaia Street, Office 54, Moskovskaia oblast, Shchyolkovo, 141109 Russian Federation. Phone +7 499 502 79 44, fax +7 499 502 50 36.
Email: info@maxmedical.ru .

If any problems with the cryosurgical instrument CryoPro occur within the warranty period please tell the model name, its serial number, the date of purchase and the description of the failure to the authorized representative of Cortex Technology ApS in Russia, MAXmedical, LLC company at the following address:

6 Sirenevaia Street, Office 54, Moskovskaia oblast, Shchyolkovo, 141109 Russian Federation.
Phone +7 499 502 79 44, +7 499 713 35 55, fax +7 499 502 50 36. Email: info@maxmedical.ru .

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