medArtis AG
Gabriel-Max-Str. 3
D-81545 München
Tel: 0049 89 6423 088
Fax: 0049 89 6423 487

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
This work is subject to copyright. All rights are reserved, whether the whole or part of the material is concerned, specifically the right of translation, reprinting, reuse of illustrations, recitation, broadcasting, reproduction on microfilm or in any other way, and storage in databanks. Duplication of this publication or parts thereof is permitted only under the provisions of the German Copyright law in its current version, and permission must always be obtained in written from medArtis AG, Munich. Violations are liable for prosecution under the German Copyright Law.
© medArtis AG, Munich 2012

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D - 81545 München/Germany

Table of Contents

Description, Composition and Structure.....	18
Indications and Applications	18
Contraindications, Precautions and Warnings	21
Operating Instructions.....	22
Side Effects	25
General Risks for Use.....	25
Special Warning	26
Interactions.....	26
Single Use, Sterilization, Shelf life and Storage	26
Forms	27
EN 980 Symbols and description.....	30

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Syntricer[®] is a ceramic bone void filler which comprises at least $\geq 95\%$ of a β -TCP (tricalciumphosphate) $\text{Ca}_3(\text{PO}_4)_2$. The composition of **Syntricer**[®] is in accordance with the ASTM F 1088 norm. **Syntricer**[®] exists as a 3-D shell scaffold that is fully interconnected.

Syntricer[®] is available in geometric shapes such as cylinders, blocks, spheres (**Syntricer**[®] Ceraball[®]), granules and rounded granules (**Syntricer**[®] Ceround[®]). **Syntricer**[®] Ceraball[®] represent micro-chambered, shell-like structured beads with a fully interconnected system of marrow spaces to form an entire implant, thus creating an osteoconductive structure for yielding a primary formation of a cancellous bone scaffold.


All geometric shapes of **Syntricer**[®] are offered with standardized micro- and macroporosities.

Syntricer[®] having a trabecular structure is offered with macropores of 300, 600 and 1000 μm . Overall, fully interconnected porosity reaches 80 – 88%.

Indications and Applications

Sintered β -TCP ceramics are applied as bone void fillers for replacing damaged, lost or removed cancellous bone, i.e. the iliac crest, in all parts of the skeleton, preferably in the metaphyseal or epiphyseal segments. **Syntricer**[®] is also used to fill bone cysts.

It is therefore used in Traumatology, Orthopaedic Surgery, Reconstructive Surgery, Maxillo-Facial Surgery and Oto-Rhino-

Laryngology, i.e. for reconstruction of the mastoid process. The material is used for all kinds of bone defects where shear, tension and bending as well as torsion and high compressive forces can be ignored. Ceramic disease, caused by relative micro-movements of the ceramic particles or by stiff implants (screws, metal implants etc.), can be avoided with a combination of materials providing damping properties, i.e. morsellized bone preferably of the same size or by a recipient bed with an intact cancellous bone environment. This results in the so-called "sandwich implantation" as a cancellous bone-chip ceramic implant: a combination of intact cancellous bone chips with ceramic granulated material. **Syntriker**[®] best fulfills its osteoconductive function in a well-vascularized recipient bed with damping properties or in combination with cancellous bone grafts.

Syntriker[®] *Ceraball*[®] with, for example, diameters of 6 mm and 4 mm, is preferably used as a bone void filler where it yields an osteoconductive structure for a physiological scaffold of cancellous bone for primary bone formation in a lamellar fashion. As a chambered material with standardized interconnections and strong capillary forces, **Syntriker**[®] *Ceraball*[®] represents the gold standard carrier for pluri-potential mesenchymal cells, adult stem cells and bone marrow, or also cells cultured by tissue engineering.

Syntriker[®] *Granules* are preferably used in combination with morsellized cancellous bone of the same size to fill bone defects or to reconstruct bone structures such as, for example, the jaw or maxillary bone ridge.

E

20

Synticer® Ceround® represents rounded granules which can be used in combination with PMMA bone cements, i.e. in a 50:50 % mixture with autologous or bank bone of the same size as described above to reconstruct, for example, the bony bed in revision surgery. The bone chips act, among other properties, as a material with damping properties (sandwich principle). A mixture of **Synticer®**, a β -TCP ceramic material, with **Synthacer®**, a hydroxyapatite ceramic material, does not influence the material properties of the single product as described and can, on the one hand, serve as an implant with damping properties (TCP) and, on the other hand, provide slow degradation of the material for long-lasting tangential adherence of bone at the cement implant (HA). A 50:50 % mixture of TCP (**Synticer®**) with HA (**Synthacer®**) is preferably used for this indication in the form of a combination of "**Synticer® Ceround®**" and "**Synthacer® Ceround®**" material.

Synticer® Cylinders are offered with different diameters and a standard length of 25 mm (15 mm for the smaller ones, respectively). The ceramic cylinder-implants represent custom-made bone void fillers for donor defects corresponding precisely to diamond-coated instruments (SDI®/Diamond TwInS™) used for joint reconstructions. The cylinders are press-fit inserted. They arrest bleeding immediately and reconstruct the cancellous framework in a physiological manner.

Synticer® Cylinders are used as a solid implant for press-fit insertion and as a bone void filler to hermetically close defects with the ceramic cover, i.e. towards the joint cavity to avoid bleeding or the loss of marrow or stem cells. A cylinder

with a closed surface on one side can be used in the same way. In both cases, the wet grinding diamond (SDI®/Diamond TwInS™) tool provides a precise bone defect for a tight press-fit closure.

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 **Syntricer®** should not be applied in the case of osteomyelitis, non-vascularized bone defects or atrophic non-unions combined with an osteosynthesis. As a non-vascularized implant, it is not protected against bacterial colonization and should not be used when focal infections are suspected.

 **Syntricer®** should not be applied in the case of defects open to the articular joint space or in the active growth plate. Non-stabilized non-unions and amphiarthroses are considered as contraindications for any ceramic material.

 The material properties of a ceramic material such as **Syntricer®**, especially its brittleness, should therefore be taken into consideration. Consequently, **Syntricer®** should never be exposed to shear, bending, and tension forces or even high compressive forces or relative motion. Micro-motion or relative movements of parts of the implants against each other or against stiff implants such as screws or against a metal stem, plate or nail cause degradation of the ceramic material. This can be avoided by combining with damping materials according to the “sandwich principle”.

The combination of ceramic implants with cancellous bone chips or morsellized bone, preferably of the same size, offers advan-

tages for healing. Ceramic implants behave in the same way in a healthy well-vascularized cancellous bone environment with damping properties. Both procedures help avoid micro- and relative motion.

Syntricer® should not be applied in the case of pre-existing calcium metabolism disorder, unless specifically intended.

Use is restricted in the case of severe metabolic diseases, uncontrollable diabetes mellitus, disorders of the calcium metabolism, treatment with steroids, immunosuppressive treatment, endocrinological bone diseases, and treatment with drugs that interfere with calcium metabolism.

Operating Instructions

Syntricer® may only be employed by, or under the supervision of, surgeons or medical doctors, or under the supervision of medical professionals with experience in the required surgical techniques and the use of biomaterials. The selection of the form of application and the surgical procedure depend on the location, type and size of the defect. In general, the medical doctor should pay attention to a well-vascularized recipient bed freed from necrotic tissue and a complete filling of the bony defect to avoid any mechanical compaction of the implant.

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Syntricer® Cylinders

For joint reconstruction, the SDI®/Diamond TwInS™ procedure comprises a grafting technique where load-bearing surfaces are reconstructed with autologous cartilage-bone cylinder grafts. These are harvested using an atraumatic wet grinding technique using diamond instruments (SDI®/Diamond TwInS™) with an internal rinsing system and where the donor bed is filled by a press-fit inserted **Syntricer®** ceramic implant as a non-load bearing bone void filler. A corresponding **Syntricer® Cylinder** is offered for each diamond tool in order to ensure a precise press-fit for filling the void.

 **Syntricer® Cylinder** implants must not be hammered; they are inserted by finger only. Any mechanical compaction has to be avoided (Fig.1).

*Fig.1 The donor beds in the patellar groove are best filled by press-fit inserted **Syntricer®** Cylinders. The ceramic cylinder must not be hammered; insertion is preferably performed with the thumb only*

It is useful to prepare defects with diamond-coated instruments (SDI®/Diamond TwInS™) which enables a press-fit seat for a ceramic cylinder.

**Synticer® Granules**

Filling the bone void with **Synticer® Granules** may be preferably combined with autologous morsellized bone of the same size or even bank bone. In the case of a bone cyst, the compact wall of the cyst should be perforated towards the free medullary cavity, providing a revascularization of the cyst cavity. Combining **Synticer® Granules** with bone marrow aspirate is an advantage. In this case, ensure that all surfaces of the implant are thoroughly saturated with bone marrow aspirate, allowing approximately 10 minutes for coagulation.

A cover should close defects in a joint, preferably via a periosteum-bone graft from the iliac crest or a solid **Synticer®** ceramic cover. Care should be taken to insert the cover at the level of the cartilage base-plate.

Synticer® Ceraball®

Synticer® Ceraball® can be useful for the reconstruction of large bone defects. The newly formed scaffold around the **Ceraball®** beads comprises lamellar bone under load. The advantage of this procedure is the early colonization of the implant which follows the surfaces of the beads, yielding a primary formation of the spongy scaffold, the remodelling process of which is reduced to a minimum. From the beginning, the newly formed bone is under load whereas the ceramic material is reabsorbed simultaneously. **Synticer® Ceraball®** represents a form-giving conglomerate as long as the defect is loosely filled without any mechanical impacting during the operation (Fig. 2).

Loose filling of the defect achieves a well-structured scaffolding; if necessary, it may be helpful to gently press a wet com-

press onto the defect and to float the beads using a syringe with physiological solution in order to achieve a perfect filling of the cavity. When using an injection device, any increase in pressure that might cause an embolism has to be avoided.

*Fig.2: The defect is loosely filled with **Syntricer**® Ceraball® comprising beads measuring 4 and 6 mm in diameter.*

A mixture of 50:50% **Syntricer**® Ceraball® and **Synthacer**® Ceraball® may be beneficial, i.e. in the case of osteoporosis, providing a slower degradation of the material on the one hand and mechanoinduction of the newly formed lamellar bone on the other hand.

Side Effects

General Risks for Use

 **Syntricer**® may only be inserted by trained surgeons. Since **Syntricer**® products are implants without an own defense mechanism, there is an imminent danger of contamination by hematogenous or intraoperative infection. As in all operations of the lower extremities, thromboembolism may be another common complication. In both these circumstances, an antibiotic prophylaxis and a thromboembolic prophylaxis, preferably using heparin-based pharmaceuticals, is strongly recommended – at least during the non-weight bearing period after the operation.

Special Warning

Syntricer® products cannot assume any load-bearing functions. Additional stabilization using osteosynthetic measures may be needed depending on the type and location of the defect to be treated.

Interactions

No interactions between **Syntricer**® products and pharmaceuticals or other medical devices have been reported to date, except the direct in-contact combination with stiff implants.

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STERILE R **Syntricer**® *Cylinders* and **Syntricer**® *Blocks* are packaged double sterile (Gamma irradiation). Implants are delivered in shockproof membrane boxes.

Syntricer® *Ceraball*® and **Syntricer**® *Granules*, as well as **Syntricer**® *Ceround*®, are packaged double sterile in small ampoules ready for application.

⊗ **Syntricer**® products are intended for single use only; re-sterilization is not allowed.

⊠ **Syntricer**® products are packaged shockproof and should only be removed from the box immediately before use; the package provides damping properties and protects the implant.

⚠ Mechanical damage can infringe on the strength of the material; therefore attention has to be paid to open the membrane-box or the ampoule only shortly before use.

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 **Syntricer®** products must not be used after the expiration date.

Forms

Syntricer® is available in a variety of forms such as cylinders, blocks, balls or beads **Syntricer® (Cerabal®)** in different sizes. The surfaces can show open pores; or one or more surfaces may be "closed" or even covered by a dense compact ceramic structure. **Syntricer®** products are available as a conglomerate of shells or a conglomerate of ceramic beads with different diameters and in form of cylinders and blocks with different lengths. The macro-pores of the structured material show an average size of 300, 600 or 1000 µm or a mixture thereof. The blocks measure 5 x 5, 10 x 10, 15 x 15, 20 x 20 and 25 x 25 mm with a height of 15 mm.

Syntricer® Cerabal® is available with different diameters. The macro-pores of the structured material show an average of 300, 600 or 1000 µm or a mixture thereof.

Syntricer® is also offered as a granulated material in natural form or rounded (**Syntricer® Ceround®**) with graduated particle sizes ranging from 0.1 – 0.9 mm (powder); 1 – 2.5 mm (medium coarse) up to 2.6 – 4.8 mm (coarse).

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This work is subject to copyright. All rights are reserved, whether the whole or part of the material is concerned, specifically the right of translation, reprinting, reuse of illustrations, recitation, broadcasting, reproduction on microfilm or in any other way, and storage in databanks. Duplication of this publication or parts thereof is permitted only under the provisions of the German Copyright law in its current version, and permission must always be obtained in written from medArtis AG, Munich. Violations are liable for prosecution under the German Copyright Law.

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88


Table of Contents

Description, Composition and Structure.....	18
Indications and Applications	18
Contraindications, Precautions and Warnings	21
Operating Instructions.....	22
Side Effects	25
General Risks for Use.....	26
Special Warning	26
Interactions.....	26
Single Use, Sterilization, Shelf Life and Storage	27
Forms	28
EN 980 Symbols and Description	29

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Description, Composition and Structure

Synthacer® is a ceramic bone void filler which comprises at least $\geq 95\%$ of a HA hydroxyapatite $\text{Ca}_5(\text{PO}_4)_3\text{OH}$.

The composition of **Synthacer®** is in accordance with the ASTM F 1185 norm. **Synthacer®** exists as a 3-D shell scaffold that is fully interconnected.

Synthacer® is available in geometric shapes such as cylinders, blocks, spheres (**Synthacer® Ceraball®**), granules and rounded granules (**Synthacer® Ceround®**). **Synthacer®** represent micro-chambered, shell-like structured beads with a fully interconnected system of marrow spaces to form an entire implant, thus creating an osteoconductive structure for yielding a primary formation of a cancellous bone scaffold.

All geometric shapes of **Synthacer®** are offered with standardized micro- and macroporosities.

Synthacer® having a trabecular structure is offered with macropores of 300, 600 and 1000 μm . Overall, fully interconnected porosity reaches 80 – 88%.

Indications and Applications

Sintered HA ceramics are applied as bone void fillers for replacing damaged, lost or removed cancellous bone, i.e. the iliac crest, in all parts of the skeleton, preferably in the metaphyseal or epiphyseal segments. **Synthacer®** is also used to fill bone cysts.

It is therefore used in Traumatology, Orthopedic Surgery, Reconstructive Surgery, Maxillo-Facial Surgery and Oto-Rhino-Laryngology, i.e. for reconstruction of the "mastoid process". **Synthacer**[®] is also used as part of an artificial orbital implant.

The material is used for all kinds of bone defects where shear, tension and bending as well as torsion and high compressive forces can be ignored. Ceramic disease, caused by relative micro-movements of the ceramic particles or by stiff implants (screws, metal implants etc.), can be avoided with a combination of materials providing damping properties, i.e. morsellized bone preferably of the same size or by a recipient bed with an intact cancellous bone environment. This results in the so-called "sandwich implantation" as a cancellous bone-chip ceramic implant: a combination of intact cancellous bone chips with ceramic granulated material.

Synthacer[®] best fulfills its osteoconductive function in a well-vascularized recipient bed with damping properties or in combination with cancellous bone grafts.

Synthacer[®] *Ceraball*[®] with, for example, diameters of 6 mm and 4 mm, is preferably used as a bone void filler where it yields an osteoconductive structure for a physiological scaffold of cancellous bone for primary bone formation in a lamellar fashion.

As a chambered material with standardized interconnections and strong capillary forces, **Synthacer**[®] *Ceraball*[®] represents the gold standard carrier for pluri-potential mesenchymal cells, adult stem cells and bone marrow, or even for cells cultured by tissue engineering.

Synthacer® Granules are preferably used in combination with morsellized cancellous bone of the same size to fill bone defects or to reconstruct bone structures such as, for example, the jaw or maxillary bone ridge.

Synthacer® Ceround® represents rounded granules which can be used in combination with PMMA bone cements, i.e. in a 50:50 % mixture with autologous or bank bone of the same size as described above to reconstruct, for example, the bony bed in revision surgery. The bone chips act, among other properties, as a material with damping properties (sandwich principle).

A mixture of **Synthacer®**, a HA ceramic, with **Syntricer®**, a β -TCP ceramic, does not influence the material properties of the single product as described and can, on the one hand, serve as an implant with damping properties (β -TCP) and, on the other hand, provide slow degradation of the material for long-lasting tangential adherence of bone at the cement implant (HA); a 50:50 % mixture of TCP (**Syntricer®**) with HA (**Synthacer®**) is preferably used for this indication in the form of a combination of "**Synthacer® Ceround®**" and "**Syntricer® Ceround®**" material.

Synthacer® Cylinders are offered with different diameters and a standard length of 25 mm (15 mm for the smaller ones, respectively). The ceramic cylinder-implants represent custom-made bone void fillers for donor defects corresponding precisely to diamond-coated instruments (SDI®/Diamond TwInS™) used for joint reconstructions. The cylinders are press-fit inserted. They arrest bleeding immediately and reconstruct the cancellous framework in a physiological manner.

Synthacer® Cylinders are used as a solid implant for press-fit insertion and as a bone void filler to hermetically close defects with the ceramic cover, i.e. towards the musculature in the iliac crest defects to avoid bleeding or loss of marrow or stem cells. A **Synthacer® Cylinder** with a closed surface on one side can be used in the same way. In both cases, the wet grinding diamond (SDI®/Diamond TwInS™) tool provides a precise bone defect for a tight press-fit closure.

Contraindications, Precautions and Warnings

 **Synthacer®** should not be applied in the case of osteomyelitis, non-vascularized bone defects or atrophic non-unions, even when combined with osteosynthesis. As a non-vascularized implant, it is not protected against bacterial colonization and should not be used when focal infections are suspected.

 **Synthacer®** should not be applied in the case of defects open to the articular joint space or in the active growth plate. Non-stabilized non-unions and amphiarthroses are considered as contraindications for any ceramic material.

 The material properties of a ceramic HA such as **Synthacer®**, especially its brittleness, should therefore be taken into consideration. Consequently, **Synthacer®** should never be exposed to shear, bending, and tension forces or even high compressive forces or relative motion. Micro-motion or relative movements of parts of the implants against each other or against stiff implants such as screws or against a metal stem, plate or nail cause degradation of the ceramic material. This can be avoided by combining with damping materials according to the “sandwich principle”.

The combination of ceramic implants with cancellous bone chips or morsellized bone, preferably of the same size, offers advantages for healing. Ceramic implants behave in the same way in a healthy well-vascularized cancellous bone environment with damping properties. Both procedures help avoid micro- and relative motion.

 **Synthacer®** should not be applied in the case of pre-existing calcium metabolism disorder, unless specifically intended. Use is restricted in the case of severe metabolic diseases, uncontrollable diabetes mellitus, disorders of the calcium metabolism, treatment with steroids, immunosuppressive treatment, endocrinological bone diseases, and treatment with drugs that interfere with calcium metabolism.

Operating Instructions

 **Synthacer®** may only be employed by, or under the supervision of, surgeons or medical doctors, or under the supervision of medical professionals with experience in the required surgical techniques and the use of biomaterials. The selection of the form of application and the surgical procedure depend on the location, type and size of the defect. In general, the medical doctor should pay attention to a well-vascularized recipient bed freed from necrotic tissue and a complete filling of the bony defect to avoid any mechanical compaction of the implant.

When filling the defect with rods or blocks, attention should be paid to the damping properties of the recipient bed or a combination with cancellous chips in order to avoid relative motion and resulting degradation of the material.

Synthacer® Cylinders

For joint reconstruction, the SDI®/Diamond TwInS™ procedure comprises a grafting technique where load-bearing surfaces are reconstructed with autologous cartilage-bone cylinder grafts. These are harvested using an atraumatic wet grinding technique using diamond instruments (SDI®/Diamond TwInS™) with an internal rinsing system and where the donor bed is filled by a press-fit inserted **Syntricer® Cylinder** implant as a non-load bearing bone void filler. As an exception, a **Synthacer® Cylinder** can be used and press-fit inserted as long as the defect is closed towards the joint using a periosteum-bone cylinder. Pronounced osteoporosis may present an indication for which **Synthacer®** has desirable properties.

A corresponding **Synthacer® Cylinder** is offered for each diamond tool in order to ensure a precise press-fit for filling the void.

 **Synthacer® Cylinder** implants must not be hammered; they are inserted by finger only. Any mechanical compaction has to be avoided (Fig.1).

*Fig.1 The donor beds in the iliac crest are best filled by press-fit inserted **Synthacer®** Cylinders. The ceramic cylinder must not be hammered; insertion is preferably performed with the thumb only.*

It is useful to prepare defects with diamond-coated instruments (SDI®/Diamond TwinS™) which enables a press-fit seat for a ceramic cylinder.

Synthacer® Granules

Filling the bone void with **Synthacer® Granules** may be preferably combined with autologous morsellized bone of the same size or even bank bone. In the case of a bone cyst, the compact wall of the cyst should be perforated towards the free medullary cavity, providing a revascularization of the cyst cavity.


Combining **Synthacer® Granule** with bone marrow aspirate is an advantage. In this case, ensure that all surfaces of the implant are thoroughly saturated with bone marrow aspirate, allowing approximately 10 minutes for coagulation.

A cover should close defects in a joint, preferably via a periosteum-bone graft from the iliac crest or a solid **Synthacer®** ceramic cover. Attention should be applied to insert the cover at the level of the cartilage base-plate. **Synthacer®** HA must not be used in donor defects of the joint except for intended reasons, i.e. in osteoporotic bone.

Synthacer® Ceraball®

Synthacer® Ceraball® can be useful for the reconstruction of large bone defects. The newly formed scaffold around the **Ceraball®** beads comprises lamellar bone under load. The advantage of this procedure is the early colonization of the implant which follows the surfaces of the beads, yielding a primary formation of

the spongy scaffold. Consequently, the remodelling process is reduced to a minimum. The newly formed bone is under load whereas the ceramic material is reabsorbed simultaneously although this is very slow as far as **Synthacer**[®] is concerned. **Synthacer**[®] **Ceraball**[®] represents a form-giving conglomerate as long as the defect is loosely filled without any mechanical impacting during the operation (Fig. 2).

*Fig. 2:
The defect is loosely filled with
Synthacer[®] **Ceraball**[®] comprising
beads measuring 4 mm and
6 mm in diameter.*

A mixture of 50:50% **Synthacer**[®] **Ceraball**[®] and **Syntricer**[®] **Ceraball**[®] may be beneficial, i.e. in the case of osteoporosis, providing a slower degradation of the material on the one hand and mechanoinduction of the newly formed lamellar bone on the other hand.

Loose filling of the defect achieves well-structured scaffolds; it may be helpful to use a wet sponge gently pressed onto the defect whereas a syringe with physiological solution is used for floating the beads to achieve a perfect filling of the cavity. When using an injection device, any increase in pressure, which might produce an embolism, has to be avoided.

Side Effects

General Risks for Use

 **Synthacer**[®] may only be inserted by trained surgeons. Since **Synthacer**[®] products are implants without an own defense mechanism, there is an imminent danger of contamination by hematogenous or intraoperative infection. As in all operations of the lower extremities, thromboembolism may be another common complication. In both these circumstances, an antibiotic prophylaxis and a thromboembolic prophylaxis, preferably using heparin-based pharmaceuticals, is strongly recommended – at least during the non-weight bearing period after the operation.

Special Warning

 **Synthacer**[®] products cannot assume any load-bearing functions. Additional stabilization using osteosynthetic measures may be needed depending on the type and location of the defect to be treated.

Interactions

No interactions between **Synthacer**[®] and pharmaceuticals or other medical devices have been reported to date, except the direct in-contact combination with stiff implants.

Single Use, Sterilization, Shelf Life and Storage

STERILE R **Synthacer®** *Cylinders*, **Synthacer®** *Blocks*, **Synthacer®** *Wedges* are packaged double sterile (Gamma irradiation). Implants are delivered in shockproof membrane boxes.

Synthacer® *Ceraball®* and **Synthacer®** *Granules*, as well as **Synthacer®** *Ceround®*, are packaged double sterile in small ampoules ready for application.

 **Synthacer®** products are intended for single use only; re-sterilization is not allowed. **Synthacer®** and should only be removed from the box immediately before use; the package provides damping properties and protects the implant

 Mechanical damage can infringe on the strength of the material; therefore attention has to be paid to open the membrane box or the ampoule only shortly before use.

 The product must not be used if the sterile packaging shows visible signs of damage. Before use, the intactness of the packaging has to be checked.

 The medical device must be stored dry in its transport box at room temperature.

 **Synthacer®** products must not be used after the expiration date.

Forms

Synthacer® is available in a variety of forms such as cylinders, wedges, blocks, balls or beads (**Synthacer® Ceraball®**) in different sizes. The surfaces can show open pores; or one or more surfaces may be “closed” or even covered by a dense compact ceramic structure. **Synthacer®** products are available as a conglomerate of shells or a conglomerate of ceramic beads with different diameters and in form of cylinders and blocks with different lengths. The macro-pores of the structured material show an average size of 300, 600 or 1000 µm or a mixture thereof.

The blocks measure 5 x 5, 10 x 10, 15 x 15, 20 x 20 and 25 x 25 mm with a height of 15 mm.

Synthacer® Ceraball® is available with different diameters. The macro-pores of the structured material show an average size of 300, 600 or 1000 µm or a mixture thereof.

Synthacer® is also offered as a granulated material in natural form (**Synthacer® Granules**) or rounded (**Synthacer® Ceround®**) with graduated particle sizes ranging from 0.1 – 0.9 mm (powder); 1 – 2.5 mm (medium coarse) up to 2.6 – 4.8 mm (coarse).

EN 980 Symbols and Description

CONSULT INSTRUCTIONS FOR USE

CAUTION

DO NOT REUSE

USE BY

BATCH CODE

CATALOG NUMBER

METHOD OF STERILIZATION USING RADIATION

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Medizinprodukte und Forschung AG
KARL STORZ GmbH & Co. KG