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
Instruction Manual

CE 0197

Store this instruction manual in a convenient location for future reference whenever necessary. Read the instructions carefully before using the product, and operate in accordance with instructions. Perform maintenance and inspections to ensure safe use and extend the life of this product.

 **TERUMO[®]**

Contents

Read this instruction manual carefully to fully understand this product and ensure its safe use. For an explanation regarding the introduction of this product, consult with TERUMO trained service technicians. Conduct safety training sessions for this product before use.

Contents

Purpose, Overview and Features of the Product —	3	Changing the Sound Volume.....	72
Purpose and Overview	3	Changing the Brightness	75
Features	3	Maintenance Timer Function	77
Unit Correlation	3	Setting the Date and Time	79
Parts Description —	4	Profile (TE-SS800 only).....	83
Exterior Diagram	4	Communication Environment (TE-SS800 only)....	85
Operation Panel	5	About External Communication Functions —	89
Screen Description	6	Overview	89
Standard Accessories	12	Network System Connection	90
Optional Accessories	12	Procedure for Connecting to Nurse Call.....	98
Precautions —	14	After Use —	99
Use	14	Cleaning	99
Compatible Syringes —	18	How to Clean Components	100
Preparation —	19	Storage.....	102
Prior to the Pump being Used	19	Maintenance and Inspections —	103
Prior-to-use Inspection	19	Maintenance and Inspection Items by TERUMO	
Attaching/Detaching a Pole Clamp.....	21	Certified Service Technicians	103
Attaching/Detaching to the IV Pole.....	23	Waste and Recycle.....	103
Connecting to AC Power	24	Troubleshooting —	104
Operating by the Internal Battery	25	Specifications —	117
Operation Procedure —	26	Units Available	122
Turning the Power On	26	Device Characteristics —	123
Setting the Syringe	27	Flow Rate Characteristics	123
Setting Flow Rate	31	Occlusion Characteristics	124
Setting VTBI	32	For Medical Staff —	125
Setting VTBI Time	33	EMC Technical Information	125
Priming	34	Symbols —	129
Inserting a Needle	35	Memo —	130
Starting Solution Delivery	36	FOR INFORMATION ABOUT TERUMO	
Stopping Solution Delivery	37	PRODUCTS —	132
Turning the Power Off	38		
Other Operation Procedures —	39		
Clearing the Volume Delivered.....	39		
Standby Function	40		
Replacing Syringe to Continue Solution Delivery ..	41		
Selecting Dose Mode (Weight Mode)			
(TE-SS800 only).....	42		
Selecting Dose Mode (Drug Library) (TE-SS800			
only).....	45		
Setting Special Functions.....	51		
Changing Settings on the Menu —	52		
Occlusion Detection Pressure	53		
Keypad Lock Function.....	56		
Bolus	59		
Advanced Dose Mode (TE-SS800 only)	65		
History Function	69		

Purpose, Overview and Features of the Product

Purpose and Overview

The TERUFUSION Syringe Pump Type SS is intended for the infusion of general and local anesthetics, vasodilators, cardiotonics, parenteral feeding, blood transfusions, anticancer drugs and labor-inducing drugs at ICU, CCU, NICU, operating rooms and general wards.

It is designed for accuracy and ease of operation, offering a precise microcomputer-controlled infusion, necessary alarms, and splash-proof construction to prevent the entry of solutions.

Model	Catalogue number	Note1	Note2
TE-SS700	TE-SS700xxx	—	Hereinafter referred to as TE-SS700
	TE-SS702xxx	With external communication (RS-232C) function	Hereinafter referred to as TE-SS702
TE-SS800	TE-SS800xxx	With infrared communication (IrDA) and wireless LAN functions	Hereinafter referred to as TE-SS800

(xxx = EN1, EN2, DE1, FR1, IT1, IT3, NL1, RU1, ES1, SV1, NO1, DA4, PL1, HU1, SK1)

Features

Dosage unit (TE-SS800 only)

- The flow rate (mL/h) can be automatically calculated based on the dosage units such as $\mu\text{g/kg/min}$. The calculated contents of the flow rate can be checked on the screen.

Easy-to-read display

- With a 4.3 inch (10.9 cm) colour LCD, clear and precise visibility has been achieved.
- The operation indicator allows the operation/stop/alarm status to be easily recognised.
- The internal battery charge and level indicator is displayed on the LCD with 5 bars.

Wide range of alarm functions

- The alarm functions include Occlusion alarm, Pressure alarm^{*1}, Nearly Empty alarm, Slider Displacement alarm, Syringe Barrel Detection alarm, Syringe Displacement alarm, Link Interruption alarm^{*1}, Plunger Displacement alarm, Battery alarm, Re-alarm, Start Reminder, No Flow Rate alarm, No VTBI alarm^{*2}, Flow Rate/VTBI (Volume to be infused) Volume Judgment alarm^{*2} ^{*3} and Completion alarm^{*2}.

^{*1}: TE-SS800 only (Default: disabled). Contact TERUMO trained service technicians for details.

^{*2}: Only if the VTBI setting function is enabled (Default: disabled)

^{*3}: The Flow Rate/VTBI Volume Judgment alarm works as a safety function to prevent the normal mode operation if the flow rate is equal to or more than the VTBI.

- The Occlusion alarm detection sensitivity level can be changed according to the condition of use.

Additional safety functions

- The keypad lock function disables any switch operations after being set to prevent any operational errors or inadvertent operations. (Except for Start switch, Power switch and Back/Mute switch (for operation to release keypad lock).)

User-friendly shape and structure

- With a unit weight of approx. 2.0 kg, the carrying burden is reduced, and allows for stable installation with the pole clamp included.
- The rounded-cornered body is easy to clean even when drug solution has adhered.

Dual power supply system

- Dual power supply system with AC power and internal battery.
- The internal battery can provide approx. 12 hours of continuous operation for solution delivery at a flow rate of 5 mL/h at ambient temperature of 25°C, with a fully-charged new battery.

Communication function (TE-SS800 only)

- The infrared communication allows communication with the optional Communication Rack System^{*4}.
- The wireless LAN allows communication with the network system.

^{*4}: TERUFUSION Communication Rack System (Model: TE-RS800)

TERUFUSION Communication Rack System (Extension) (Model: TE-RS811)


Unit Correlation

mcg/kg/min and $\mu\text{g/kg/min}$ have the same meaning and are interchangeable.

Parts Description

Exterior Diagram


◆ Front view


◆ Rear view

TE-SS700
TE-SS800

TE-SS702

Operation Panel

Name	Function
 Power switch	Turns the power ON/OFF.
 Purge switch	In the stop status, press and hold this switch to perform a rapid infusion for priming. In the bolus (Hands On/Hands Free) mode, press or press and hold the switch to start the bolus infusion. The buzzer sounds and the operation indicator flashes in green.
 Start switch	Starts the solution delivery.
 Stop switch	Stops the solution delivery. When setting the flow rate etc., rotating the Dial while pressing the Stop switch allows changing the setting at the rate of 100 times.
 Display select switch	Switches the selected item (Flow rate screen, Menu screen). Press for 2 seconds or more on the flow rate screen to go to the menu screen.
 Back/Mute switch	Mute Cancels the confirmation. Returns to the previous menu screen (the setting will not be cancelled). In the stop status, press and hold for 2 seconds or more to standby. In the keypad lock status, press and hold for 2 seconds or more to release the keypad lock.
 Confirmation switch	Confirms the item selection and is also used for check confirmation.
 Operation indicator	Displays the operation status. Green flashing: Solution delivery in progress No light: Stopped Red flashing: Stopped by an alarm Green and red flashing alternately: Self-check (self-diagnosis) operation in progress / Standby mode Red light: Device fault


Parts Description

Screen Description

Note

- As an example, this manual uses screen display examples from TE-SS800.


◆ Flow rate screen


1	Flow rate display	Displays the value and information for the flow rate (mL/h), etc.
2	Volume delivered display	Displays the volume delivered.

Notes

- Each time you press the Display select switch, active items rotate as shown below. Active item turns dark blue.
- Returns to the flow rate screen after 10 seconds of no operation.

Status icon display

 Battery icon Displays the remaining battery level.	 (lit)	When AC power is connected and charge is full.
	Displays  and the charge level alternately.	When AC power is connected and battery is charging.
	 (flashing)	When the battery is not connected or faulty. (Replacement cycle of the battery is 2.5 to 3 years)
	Battery charge level (displayed in 5 bars)  5 bars (green): Approx. 12 hours  4 bars (green): Approx. 8 hours  3 bars (green): Approx. 6 hours  2 bars (green): Approx. 4 hours  1 bar (red): Approx. 30 mins (when the Battery alarm has not been issued) (In the case of a continuous solution delivery at a flow rate of 5 mL/h at ambient temperature of 25°C, with new battery after charging for 8 hours or more with the power turned off.)	
 AC icon Displays connection/disconnection of the AC power supply.		When AC power is not connected
		When AC power is connected
 Occlusion icon The occlusion detection pressure is displayed in 10 levels (1 to 10) in the upper right of the icon.		The occlusion detection pressure level is 1.
		The occlusion detection pressure level is 6 (Default).
		The occlusion detection pressure level is 10.
 Occlusion icon The internal pressure level of the infusion line that does not reach the level of Occlusion alarm is displayed in 4 levels.	 (grey)	Below 25% of the specified occlusion detection pressure.
	 (green)	25% or more of the specified occlusion detection pressure.
	 (yellow)	50% or more of the specified occlusion detection pressure.
	 (orange)	75% or more of the specified occlusion detection pressure.

Parts Description

 Syringe brand/Syringe size icon		Syringe brand	
		Syringe size to be used (5, 10, 20, 30, 50 mL) Note • If a syringe has not been set or has not been correctly detected, the size display indicates 	
 Keypad lock icon Disables the switches to prevent an inadvertent operation.		When keypad lock is enabled (ON).	
	(No indicator)	When keypad lock is not disabled (OFF).	
 Communication icon Displays during infrared communication (IrDA) or external communication (RS-232C).		During communication	
 Wireless LAN icon Displays the strength of the signal. (TE-SS800 only)		Antenna only is flashing: The wireless LAN module is faulty. Antenna only is lit: Connection to the network system is not established.	Antenna + signal strength bars light grey: Connection to the network system is established. Antenna + signal strength bars light green: During communication.
		Signal condition: Poor	
			
			
			
			
 Maintenance icon Displays when the maintenance timer set period has elapsed.		When the set time period has elapsed	

◆ Sensor icon display

The icons except for the occlusion icon will flash when the syringe is disconnected regardless of an alarm event.

Syringe Displacement icon	Displays when there is a Syringe Displacement alarm.
Nearly Empty icon	Displays when there is a Nearly Empty alarm or Syringe Empty.
Occlusion icon	Displays when there is an Occlusion alarm or Syringe Empty.
Slider Displacement icon	Displays when there is a Plunger Displacement alarm or a Slider Displacement alarm.

Note

- Nearly Empty icon and Occlusion icon are displayed when Syringe Empty occurs. (See page 109 "Syringe Empty")

◆ Tab display

* A selected item that can be entered is displayed in dark blue.

 Flow rate tab	Used for setting the flow rate.
 Flow rate tab (TE-SS800 only)	Used for setting the flow rate when the library or the weight mode is selected.
 Volume delivered clear tab	Used for clearing volume delivery.
 Info tab	Displays the setting information (unit, weight, dilution and etc.).
 Hands On Bolus tab	Used for setting the Hands On Bolus function.
 Hands Free Bolus tab	Used for setting the Hands Free Bolus function.

◆ Ready icon


 Start up icon	Displays when solution delivery is ready to start.
 Ready for bolus icon	Displays when the pump is ready to start bolus.

Parts Description

◆ Menu screen

For the programming procedure, see page 52 "Changing Settings on the Menu".


Menu 1

Menu tab


Menu details display


Menu 2

Menu tab

Menu details display

Menu	Menu tab	Menu details	Note
Menu 1 Display select switch (Press and hold for 2 seconds or more) 	 Dose mode select (TE-SS800 only)	mL/h µg/kg/min mg/kg/h Library ("Library" is displayed only after uploading the library.)	Other dose units can be specified in the Drug Library.
	 Occlusion detection pressure	1 to 10 (± 10 kPa ↔ ± 120 kPa)	Sets the occlusion detection pressure level.
	 Keypad lock	ON OFF	Sets the keypad lock function to ON or OFF.
	 Bolus	Hands On Bolus Hands Free Bolus	Bolus infusion is available with the Purge switch.
	 Advanced dose mode (TE-SS800 only)	OFF Interval/Intermittent Mode Multi-step Mode Delayed Start Mode	Selects an advanced dose mode. Displayed when a drug library that allows for advanced dose mode is selected.
	 Combination dose mode (TE-SS800 only)	—	For details, contact TERUMO trained service technicians.


Menu	Menu tab	Menu details	Note
Menu 2 Press and hold Back/Mute switch, then press and hold Display select switch for 2 seconds or more. 	 History	up to 10,000 events	Checks the operation history.
	 Sound volume	Alarm sound: 1: Low 2: Medium 3: High Step: Step up Operation sound: 0: Very low 1: Low 2: Medium 3: High Cable insertion or extraction sound: 0: Very low 1: Low 2: Medium 3: High	Sets the sound volume level.
	 Brightness	1 to 5 (Dark ↔ Bright)	Sets the brightness level.
	 Maintenance timer	OFF to 36 months	Sets a period before the maintenance timer display is lit.
	 Date and time	Year, Month, Day, Hour, Minute	Sets the date and time.
	 Profile (TE-SS800 only)	Profile range can be set through TERUFUSION Drug Library Manager or TERUFUSION Software Package.	Selects the profile Up to 30 items
	 Communication environment (TE-SS800 only)	Selection of Communication environment table (up to 10) Wireless LAN ON/OFF	Sets the communication environment.
	 Password	If the password is set, some special functions are password protected.	Displays the menu only available to TERUMO trained service technicians.

Parts Description


Standard Accessories

Standard accessories supplied with this product:


- AC power cable

- Pole clamp (Screw type)

- Instruction manual


Optional Accessories

Optional accessories are sold separately.


The specifications and external appearance of the product may be changed without notice for the purpose of improvement.

Name	Model	Catalogue number
TERUFUSION One Touch Pole Clamp	TE-877	TE*877
TERUFUSION Drug Library Manager	TE-SW800B	TE*SW800BE
TERUFUSION Software Package	TE-SW800P	TE*SW800PE
TERUFUSION Standard Rack System	TE-RS700	TE*RS700N
TERUFUSION Communication Rack System	TE-RS800	TE*RS800N TE*RS800N03
TERUFUSION Communication Rack System (Extension)	TE-RS811	TE*RS811N


- TERUFUSION One Touch Pole Clamp (TE-877)


- TERUFUSION Drug Library Manager (TE-SW800B)
(TE-SS800 only)
See the instruction manual for details.


- TERUFUSION Software Package (TE-SW800P)
(TE-SS702/TE-SS800 only)
See the instruction manual for details.


- Rack system
TERUFUSION Standard Rack System (TE-RS700)


- TERUFUSION Communication Rack System (TE-RS800)

<Combination Examples>

- TERUFUSION Communication Rack System (Extension) (TE-RS811)

Notes

- See each rack system instruction manual for details.
- Up to three rack systems can be joined for use.
- To join rack systems, contact TERUMO trained service technicians.

Precautions

To ensure safe and correct use of this product, please observe all precautions.

Non-compliance with precautions and incorrect use may result in damage or injury.

The following are signs used in this manual and their meanings:

Warning

This label preceding a precaution indicates that there will be a possible risk of death or personal injury if the precaution given is not complied with.


Caution

This label preceding a precaution indicates that there will be a possible risk of personal injury or property damage if the precaution given is not complied with.

Use

Warnings

- After the power is turned on, a syringe that corresponds to the syringe brand name displayed on the LCD should be used.
[The flow rate accuracy and alarm function cannot be guaranteed if a wrong syringe is used.]
- A syringe other than specified should not be used for this product.
[If using a syringe other than those specified, the flow rate accuracy and alarm functions cannot be guaranteed.]
- When setting a syringe, ensure that the syringe plunger has been firmly set onto the slider hook and that the syringe has been correctly set into the sections including checking that the flange of the syringe has been inserted into the flange holder. The elevation difference between this product and the patient should be kept to a minimum.
[Drug solution may not be delivered correctly due to rapid infusion caused by displacement of the syringe plunger from the slider hook, incorrect detection of the syringe size, etc.]
- Always check the infusion status (decreasing volume of drug solution), connecting site and puncture site at the start of the infusion. In addition, conduct the same check during the infusion on a regular basis, including patient rounds.
- When an occlusion occurs due to tube bend in the infusion line, blocked filter, blood clot inside the needle or other reasons, close the infusion line on the downstream side.
[1. The internal pressure of infusion line from the pump to the occluded section in the downstream is high. Simply removing cause of the occlusion will result in "bolus infusion (temporally excessive infusion of drug solution)" to the patient. 2. The pump will not operate until the occlusion is removed.]
- Do not apply any pressure to the syringe plunger or slider.
[The syringe plunger may be pushed, causing bolus infusion to the patient.]
- This product should be placed and used on a level and stable location. In addition, when using the IV pole, ensure that the stand is stable and this product is securely fixed to it.
[Dropping or falling may result in damage or failure.]
- Since this product is precision equipment, it should not be used if it has received any impact (drop to floor, falling from the IV pole, violent shock).
[Even though no fault is observed in the product appearance, the original functionality or performance (flow rate accuracy and various alarm functions, etc.) of this product may not be achieved, and therefore inspection is required.]
- When using devices (mobile phones, radio devices, radio knives, defibrillators, etc.) that emit an electromagnetic wave within the area of this product, they should be used as far away as possible and you should verify the normal operation of this product in the configuration it is used. In addition, this product should use a power supply of a system isolated from these devices and should obtain secure earth connection.
[Malfunction in this product may result in critical harm to the patient.]
- No part of this product, including the infusion line, should be taken into or used in a control area for radiation devices/ MRI or inside a hyperbaric oxygen therapy room. If any part of this product is accidentally brought into these environments, immediately discontinue use.
[This product is not designed to be used in these environments. Malfunction, damage or degradation of this product may occur or it may lead to explosion.]


Warnings

- Since this product does not have an airtight structure, it should not be used or stored in an active gas environment (including sterilizer gas), nebulizer-sprayed environment, high-humidity environment, etc. It should not be submerged into water.
[If the electronic components inside the product are affected, there may be subsequent damage and time degradation which will cause failure of this product.]
- The product should not be used or stored in a flammable environment.
- This product should not be used in conjunction with gravity infusion.
[1. If this product suffers a downstream occlusion lower than the connection with the gravity infusion line, the Occlusion alarm will not operate. 2. If bubbles are entrained in the downstream of the infusion line due to emptied gravity infusion line first, the normal infusion cannot be carried out.]
- The product should not be used for an extracorporeal circuit etc. that may lead to an extreme negative or positive pressure.
[Solution delivery may not be carried out correctly or rapid infusion may occur due to displacement of the syringe gasket from the plunger or displacement of the syringe plunger from the slider hook. In addition, solution delivery may not be carried out correctly because of bolus infusion, back-flow, etc.]

Cautions

- This product should only be operated by skilled personnel.
- Before use, check the manuals of the medical supplies and medical equipment that are used in combination.
- Regular inspections should be conducted. If any fault is observed, stop using this product and request an inspection and repair. [The original functionality or performance of this product may not be achieved.]
- Prior to first use, or if unused for a long period, connect this product to an earthed AC power supply and provide a sufficient charge (8 hours or more) with the power turned off. [If not sufficiently charged, this product may not be able to operate using the internal battery during a power failure, etc.]
- Use the AC power cable included and connect to an earthed AC power source. The AC power cable included should not be used for other equipment. [Using an AC power cable other than the one specified may result in failure of this product. In addition, if used without an earth connection, the electrical safety of this product is not guaranteed.]
- Use an earthed AC power supply for normal use. The internal battery is an auxiliary power source for when the AC power supply is not adequately available during transportation, power failure, etc.
- This product should only be used with a sufficient power supply. [If the power supply is not sufficient, the internal battery will be used for the operation, which may result in no power in an emergency situation.]
- Since drug solution may cause a short circuit, ensure that the connecting sections of the AC inlet and AC power cable are not wet when connecting. If moisture is present, ensure that the power is turned off and the AC power cable is removed from both the pump and earthed AC power source, and then thoroughly wipe with a dry cloth. [Since this product does not have a waterproof structure, it may affect the electrical components inside, and cause a malfunction.]
- Attention is required not to apply any strong static electricity. [It may result in failure or malfunction.]
- When used in the vicinity of a radio knife (Medical radio knives are surgical equipment for incision and coagulation by high energy radiofrequency current), check the following before use.
 - (1) Radio knives have a different level of high-frequency noise emission depending on the types, and using old models (vacuum tube gap type) in combination should be particularly avoided as the noise levels from those are higher.
 - (2) The distance from the radio knife cord (knife holder, knife cord and return electrode cord) and radio knife body to this product should be kept to a minimum of 25 cm.
 - (3) The radio knife and this product should be operated by a power supply from a separate system, and both should be securely earthed.
- When fixing this product to the IV pole, the specified pole clamp should be used. [If a pole clamp that has not been specified is used, the function cannot be guaranteed. It may cause failure or accident.]
- When attaching a pole clamp to the pump, ensure that it is set securely. [Failure to do so may result in falling of the pump or pole clamp.]
- When a pole clamp is used, do not release it until you have checked that this product is securely fixed to the IV pole.

Cautions

- When transporting with the pump fixed to the IV pole, do not hold the handle of the pump nor apply pressure from the top. [The pole clamp may be dislocated or damaged.]
- When using the IV pole, follow the instructions for the IV pole.
- Turn the power on before the syringe is inserted, and check that the LCD and operation indicator are flashing, and that the buzzer sounds. [If the power is turned on when a syringe is installed, this product cannot operate a normal self-check (self-diagnosis).]
- When using a syringe (sterilized specified model for general use) for this product, be careful that the nominal capacity of the syringe is not exceeded. [If exceeded, the accuracy of flow rate or the alarm function cannot be guaranteed.]
- When installing a syringe, pull forward the clamp and turn. Then, after the syringe is set, turn the clamp back to the original position and gently push it to fix the syringe. In addition, the slider should be moved while holding the clutch. [Any unreasonable operation or excess operation will cause failure of this product.]
- When setting the plunger to the slider, ensure that there is no gap between the syringe plunger and slider. [If the plunger is set and there is a gap, the air or drug solution may be sucked by the slider hook.]
- After a syringe has been installed, do not apply any force such as pulling or pushing the infusion line. [Applying these forces may result in displacing the syringe barrel from the predefined position, temporarily delivering or sucking the drug solution.]
- Before inserting an intravenous needle or connecting an infusion line with the connecting site, always prime the infusion line by pressing the Purge switch. This product may cause harm to a patient due to air injection, since it does not have a function to detect air in the infusion line.
- When removing the syringe installed in this product, the three way stopcock (if used), of the infusion line should be closed before removing it. [It may result in excessive infusion (siphoning free flow by gravity) of drug solution.]
- Before starting solution delivery, check that the settings for any errors. (e.g. an error in the digit of flow rate etc.) [This product does not have a function to determine the correct values, possibly resulting in excessive infusion or insufficient infusion to a patient.]
- When starting solution delivery, check the volume delivered and then clear it as needed before use.
- When purging, check the purge volume as rapid infusion continues after the Purge switch is pressed. [It may result in excessive infusion.]
- When purging a drug solution with high viscosity through a thin intravenous needle, an Occlusion alarm is likely to be issued even if the infusion line is not occluded. A solution with high viscosity, deliver it at the rate of 100 mL/h or lower, not by purging it. [If purging is continued, an Occlusion alarm may frequently be issued or the solution may be unable to be delivered.]
- After an Occlusion alarm is issued, always eliminate the cause of occlusion before starting. [In order to correct the volume of bolus at the time of handling occlusions, this product has functions to automatically reduce the pressure inside the infusion line at the occurrence of the Occlusion alarm (pulling back the slider, subtraction from the volume delivered). Therefore, if starting without removing the cause, the infusion may not be carried out correctly such as causing the drug solution in the infusion line to flow back to the syringe, or a repetitive Occlusion alarm status.]
- When using at a low flow rate or under a low temperature, monitor the solution delivery carefully to check for any occurrence of occlusion. [The solution delivery may be interrupted for a period due to: 1. The flow rate setting decreases, the time from occlusion occurrence to detection becomes longer. 2. The temperature decreases, it interferes with the movement of the syringe (the sliding resistance of the plunger increases), frequently causing the Occlusion alarm to be issued.]
- Use the syringe after it has adjusted to the room temperature. [If the syringe is used when it is still cold, it interferes with the movement of the syringe (causing resistance of the plunger), frequently causing the Occlusion alarm to be issued.]
- Even within the normal conditions of use, this product should not be used under conditions that cause a sudden temperature change. [Condensation inside the product results in damage and time degradation, and therefore the original functionality or performance of this product may not be achieved.]
- In a case where a manufacturer has made changes to their syringe, the flow rate accuracy and alarm functions may not be guaranteed. If any fault is observed, immediately stop using the product and contact TERUMO trained service technicians.
- When connecting this product with other medical equipment or network system, check its conformance to IEC 60601-1-1: 2000 and IEC 60601-1-2: 2001 (Amd.1: 2004) prior to use in order to ensure the system safety.
- When connecting this product with a network system or other equipment, for safety, checking the manufacturer's specifications of the equipment is recommended. In addition, use EMI compliant products for the connection cable.
- When using the external communication function, pay particular attention, as it is more susceptible to the effect from a radio knife, mobile phone, radio device, defibrillator, etc. Regularly check that this product is operating normally.

Cautions

- If connecting with a nurse call, always check in advance with the electrical contractor for compatibility, and connect with a system for which the connection function has been confirmed. [The functions regarding connection are not guaranteed.] (TE-SS702 only)
- If connecting with a buzzer or lamp instead of a nurse call, it should be used with a capacity of 12V DC, 1A or less. (TE-SS702 only)
- When using the wireless LAN communication function, the effect on other equipment by radio interference should be considered. [It may affect other equipment.] (TE-SS800 only)
- When connecting this product with a network, the settings of this product and the network system need to be compatible. Contact a vendor specialist for the correct settings and allow only the system administrator to perform settings. [Failure to connect with correct settings may result in interfering with the original functionality or performance of this product and affecting the network system.] (TE-SS800 only)
- When this product is transported while in use, do not touch the switches, etc. Use the keypad lock function as needed. [Touching the switches etc. without keypad lock enabled may result in an operation (power on/off, stop, start, rapid infusion) that is not intended by this product.]
- The LCD and the Operation panel (switches etc.) of this product should not be pressed with excessive force and should not be operated with a pointed object. [It may cause damage to or failure of the LCD or the Operation panel.]
- Do not disassemble, make alterations (including actions that interfere with the functionality or performance such as taping the LCD or a movable part) to, or repair this product. [This may result in failure, damage or device performance degradation of this product.]
- The AC power cable, communication or nurse call cable, used for this product should not be pinched with forceps or equipment, e.g. caster, or punctured with a needle.
- This product should not be used in a place where vibration, dust, mist, corrosive gas, etc., occurs or in a place where the product is sprayed with liquid. If solution is spilled on this product, thoroughly wipe with a dry soft cloth.

Compatible Syringes

Use a syringe specified below.

This pump accepts TERUMO syringes and other major syringe brands listed in the table below.

Selection of a brand can be made by simple calibration. Contact TERUMO trained service technicians.

Warnings

- Before setting a syringe, ensure the brand of the syringe to be used is the same as the one specified on the pump. If the brands do not match, the flow rate accuracy and alarm functions can not be guaranteed.
- Before using the syringe pump, ensure that the size dimension of the syringe is the same as listed below. If there is any change in the syringe specification, the flow rate accuracy and alarm functions may not be guaranteed. If any fault is observed, immediately stop using the product and contact TERUMO trained service technicians.
- Only use specified syringes with luer lock for the pump.

◆ List of Compatible Syringe Brands

For your reference, these are the dimensions as measured in March 2012.

Display	Syringe Brand	Acceptable syringe sizes									
		5 mL		10 mL		20 mL		30 mL		50/60 mL	
		a [mm]	b [mm]	a [mm]	b [mm]	a [mm]	b [mm]	a [mm]	b [mm]	a [mm]	b [mm]
TERUMO_J	TERUMO JAPAN	15	14	17	18	22	19	26	23	—	—
TERUMO_P*1	TERUMO PHILIPPINES	15	14	17	18	22	19	26	23	32	24
NIPRO	NIPRO	15	14	17	18	22	19	26	23	32	24
ZENECA	ZENECA DIPRIVAN	—	—	—	—	—	—	—	—	30	17
BD_USA	BD PLASTIPAK USA	14	14	17	13	22	16	24	17	30	17
B.BR_OMN	B BRAUN OMNIFIX	14	17	18	16	22	16	24	17	31	20
B.BR_PER	B BRAUN ORIGINAL-PERFUSOR	—	—	—	—	21	16	—	—	31	38
P_FERTE	PENTAFERTE	15	16	17	17	22	16	26	13	29	18
CHIRANA	CHIRANA	—	—	—	—	—	—	—	—	30	21
CODAN	CODAN	14	14	17	13	23	17	—	—	30	20
CODAN_30	CODAN	—	—	—	—	—	—	23	17	—	—
CODAN_PR	CODAN PERFUSION	—	—	—	—	—	—	—	—	30	37
BD_PRECS	BD PRECISE	—	—	—	—	23	21	—	—	32	24
BD_PERFS	BD PLASTIPAK PERFUSION	—	—	—	—	—	—	—	—	31	38
BD_IRELD	BD PLASTIPAK IRELAND	—	—	—	—	21	16	24	16	30	16

*1 50/60 mL: The syringes made by TERUMO EUROPE can be used in this setting.

Size a : Approximate outer diameter of syringe

Size b : Approximate length from syringe flange to the end of the plunger when the plunger is completely slid into the syringe.

Cautions

- TERUMO cannot accept any responsibility for errors in flow due to modifications of the syringe specifications introduced by the manufacturer.
- Differences in factors such as size and plunger force in compatible syringes may not fulfill flow rate accuracy of infusion.
- Contact the syringe supplier for details on syringe specifications.

Note

- When a custom syringe list is created by specifying a brand for each syringe size among the available syringes, it can be registered to the pump to be available for use. Contact TERUMO trained service technicians.

Prior to the initial use of the pump, set the date and time to record history correctly. (see page 79 "Setting the Date and Time".) Before use, carefully read "Precautions" and the labels attached to the product.

Prior to the Pump being Used


Connect to the AC power supply and charge (8 hours or more) with the power turned off.

Prior-to-use Inspection

Check the following points before use.

- 1) **There is no damage to the main unit, AC power cable or pole clamp, and that there is no malfunction due to adhered drug solution.**

◆ Main unit/AC power cable


1. There is no damage, dent or adhered drug solution on the AC inlet and AC power cable.
2. There is no damage to the handle, rubber protector, tube holder and RS-232C connector cover (if RS-232C is connected) (TE-SS702 only).
3. There is no damage or adhered drug solution on the slider, slider hook and flange holder (including Slit).
4. When moving the slider left and right with the clutch pushed, it should move smoothly, and when it is moved to the left most position while the power is on, the Nearly Empty alarm should be issued.
See page 109 "Nearly Empty alarm" in "Troubleshooting".
5. There is no adhered drug solution on the syringe clamp and that it moves back and forth smoothly.
6. There is no drug solution on the Dial, it turns smoothly, and flow rate setting can be set when the power is on.
If there is adhered drug solution, promptly clean in accordance with the directions on page 100.


Preparation

Notes

- If the tube holder comes off, follow the procedure below to attach it.
- When inserting the notch of the tube holder into the hole on the rear side, insert while pulling the tube holder towards the rear side. Forcibly inserting into the hole may cause damage to the notch.

◆ Pole clamp

1. There is no damage to any part of the pole clamp and that the clamp is clean.
2. When attaching pole clamp to the pump, it should lock in place firmly with an audible click.
3. Clean periodically even if there is no adhered drug solution.

2) The following self-check (self-diagnosis) operations should be performed when the power is turned on using the internal battery without a syringe installed.

1. The operation indicator should repeatedly flash red and green.
2. The message "Self-check" is displayed on the LCD, and the LCD flashes.
3. The buzzer sounds.
4. The alarm icons, slider displacement icon , and syringe displacement icon  are flashing.

3) Battery alarm has not been issued.

1. The level of the battery icon on the LCD has two bars (green)  or more.
2. Battery alarm (yellow) is not displayed on the LCD.

4) When the AC power cable is connected, the LCD displays AC icon .

1. When the AC power cable is inserted into the AC inlet on the back of the main unit completely, the AC icon  is displayed.
2. If the internal battery is not fully charged, fully charge it.

Note

- During charging,  and  flash alternately.

5) When the AC power cable is connected, there is no looseness or any fault.


1. The AC power cable is completely inserted into the AC inlet of the main unit.

Attaching/Detaching a Pole Clamp

◆ Attach the pole clamp

The pole clamp can be locked at the pole clamp lock portion only in the direction shown in the figures below.

To attach to a vertical pole

- **Insert the pole clamp along the grooves on the back of the pump from the lower side.**


Check

- The pole clamp must be securely inserted. The pole clamp lock portion should make a click when the pole clamp is attached.

Note


- For instructions on attaching the IV pole and other parts, see page 23 "Attaching/Detaching to the IV Pole".

To attach to a horizontal pole


Notes

- The pole clamp can be attached in the vertical direction as shown in the figure. Note that the pump cannot be installed in the rack system when the pole clamp is attached in the vertical direction.
- A pump with an RS-232C connector cover (TE-SS702 only) may produce difficulties attaching to some horizontal bars depending on the diameter of the bar. (It can be installed on to "Horizontal Rail")
- Do not attach the pump in the vertical direction.


◆ Detach the pole clamp

- Pull the pole clamp downward while pressing and holding the pole clamp lock portion (green lever at the bottom of the pole clamp).

Attaching/Detaching to the IV Pole

◆ Pole clamp (Screw type)


- Turn the pole clamp knob to fix the pump firmly to the IV pole.

Note

- When removing the pump from the IV Pole, hold the main unit and turn the pole clamp knob opposite to the arrow shown in the figure.

◆ TERUFUSION One Touch Pole Clamp (Optional accessory)

- Turn the pole clamp knob to fix the pump firmly to the IV pole.

Notes

- To attach the pole clamp quickly, before fixing it to the IV pole completely, pull and hold the lever (1), push the knob in the direction of (2) and turn the knob in the direction of (3) to fix the pump. (The lever cannot be pulled after the pump is fixed to the pole).
- To detach the pole clamp from the IV pole quickly, loosen it slightly by turning the knob in the reverse direction of (3), pull and hold the lever (1), and pull the knob in the reverse direction of (2).

Warning


- Make sure that pole clamp is firmly fixed to the IV pole and the IV pole is stable. When transporting with the unit fixed to the IV pole, be aware as the pump may drop or fall due to a floor step, slope, contact with surrounding objects, etc.

Cautions


- When a pole clamp is used, do not release it until you have checked that this product is securely fixed to the IV pole.
- Do not turn the knob while the lever has been pulled. It may cause damage.
- Do not pull the lever while it has been fixed to the IV pole. It may cause damage.

Note


- The syringe pump can be set tilted with the pole clamp attached as shown below.

Connecting to AC Power

- 1) Firmly connect the AC power cable (included) to the AC inlet.
- 2) Connect the plug to an AC power source with an earth connection.

Check

- The AC icon and battery icon should be displayed on the LCD.

Notes

- Even if the power is off, if the AC power has been supplied, the AC icon is displayed.
- The battery can be charged even if the power is off. The battery icon (during recharging)  is displayed.
- If the AC power cable is disconnected or connected with power turned on, the buzzer sounds.
- If the AC power cable is disconnected with power turned on, the LCD becomes dim (2 levels of brightness setting). When the AC power cable is reconnected, the previous brightness is restored.

Cautions


- Please ensure AC inlet and all connections are clean and dry prior to use. If moisture is present, make sure that the power is turned off and the AC power cable is removed from both main unit and AC power source, and then thoroughly wipe with a dry cloth.
- Use the AC power cable included and connect to an earthed AC power source.
- Use an earthed AC power supply for normal use. The internal battery is an auxiliary power source for when the AC power supply is not adequately available during transportation, power failure, etc.
- This product should only be used with a sufficient power supply. If the power supply is not sufficient, the internal battery will be used for the operation, resulting in no power in an emergency situation.

Operating by the Internal Battery

If the AC power supply is cut off, while in use or the power supply is reduced, the power supply is automatically switched to the internal battery. The operation with internal battery is available for approx. 12 continuous hours. (In the case of continuous solution delivery at a flow rate of 5 mL/h at ambient temperature of 25°C, with new battery after charging for 8 hours or more with the power turned off.)

When operating on the internal battery, the brightness of LCD backlight becomes two levels lower.

* In the case of no AC power, the LCD displays the AC icon 

Battery charge level (displayed in 5 bars)

-  5 bars (green): Approx. 12 hours
-  4 bars (green): Approx. 8 hours
-  3 bars (green): Approx. 6 hours
-  2 bars (green): Approx. 4 hours
-  1 bar (red): Approx. 30 mins (when the Battery alarm has not been issued)

(In the case of continuous solution delivery at a flow rate of 5 mL/h at ambient temperature of 25°C, with new battery after charging for 8 hours or more with the power turned off.)

Cautions

- Prior to first use after purchase, or if unused for a long period, connect this product to an AC power supply with an earth connection and provide a sufficient charge (8 hours or more) while it is turned off. If not sufficiently charged, this product may not be able to operate using the internal battery during a power failure, etc.
- Time degradation of the internal battery will cause shorter operating hours indicated by the battery icon on the LCD. The battery should be replaced regularly approx. every 2.5 to 3 years. (Since the flashing cross on the battery icon indicates battery failure or replacement timing, use it as an indication for battery replacement.)
- If any fault is observed while operating using the internal battery, immediately stop using the product and contact TERUMO trained service technicians.

Notes

- The operating hours of the internal battery are only an estimation. Since the actual operating hours may be shorter pay extra attention to the status of the main unit when the internal battery is used.
- When the battery charge level goes down, the Battery alarm sounds.
Connect the AC cable immediately.
The alarm will continue to sound for approx. 30 minutes before the infusion is automatically stopped.
(Even after muting, the alarm will sound again every 2 minutes)
If no action is taken, a Shutdown Notice will appear and the solution delivery will be stopped. Then approx. 3 minutes later, the power is automatically turned off.
(The Shutdown Notice will continue to sound until the AC cable is connected. It cannot be muted by the Back/Mute switch.)


Operation Procedure

Turning the Power On

Please ensure the product is fit for purpose prior to use.

Cautions

- Conduct the self-check (self-diagnosis) of this product without a syringe installed and with the clamp back to the original position. Otherwise the self-check may not be conducted properly.
- Check for any fault by conducting a self-check. If a fault is observed, immediately discontinue use and contact TERUMO trained service technicians.

- Press and hold the Power switch (for 1 second or more) until the TERUMO logo screen is displayed on the LCD to turn the power on.

Startup


When turning the power on, the following occur simultaneously;

- TERUMO logo appears on the display
- LCD flashes 3 times
- Operation indicator repeats flashing red and green alternately, and then buzzer sounds.

2012/05/01 01:21

TERUMO

TERUFUSION Pump TE-SS800N
Self-check

Operation indicator

2012/05/01 01:21

TERUMO


Periodical Service Required

Notes

- When the period set by the maintenance timer is reached, a maintenance request is displayed at start up. Please contact the person responsible for servicing this product.
- After the period specified for the maintenance timer has passed, this screen displays at every start up until the setting is changed.

Checks

- The syringe displacement icon and the slider displacement icon displays and flashes on the LCD.
- The AC icon displays.
- The battery icon displays.
- The operation indicator is off.


Operation indicator

Setting the Syringe


- 1) Attach the IV line to a syringe filled with solution using aseptic technique.

Warning


- Only use specified syringes with luer lock for the pump.

- 2) Pull the clamp forward and turn.


Place the syringe


- 3) Insert the flange into the slit.

Checks

- Check the LCD as shown below.
- The syringe flange is correctly fitted into the flange holder.

Place the syringe


- 4) The syringe plunger must be held securely between the slider and the slider hook.

While holding the clutch, move the slider until it meets the plunger, and then release the clutch. If the following guidance in the figure below does not display, perform this step again.

Check

- The clutch displays in the LCD as shown below.
- The slider displacement icon is turned off.

Place the syringe



Note

- If the process is not completed after several attempts, the clutch may be ill fitting. Ensure that the IV line is not attached to a patient and then apply light pressure to the slider to push the clutch into position. The clutch may click when fitted into the position.


- 5) Turn the clamp back to the original position and gently push it to fix the syringe.

Check

- The message "Syringe Barrel OK xx mL" is displayed on the LCD, and the syringe displacement icon is not displayed.

Syringe Barrel OK 50mL



- 6) Turn the Dial or press the Display select switch to select a syringe brand.

Warning

- After the syringe type is selected, recheck that the setting has been adjusted correctly.

Notes


- Turning the Dial or pressing the Display select switch moves the white frame.
- The pump will automatically detect the syringe size.


- 7) Press the Confirmation switch.

Check

- The syringe brand and the syringe size are displayed on the LCD correctly.

- 8) Insert the tube completely into the tube guide and set the tube to the tube holder.

Operation Procedure

Warning


- When setting a syringe, ensure that the syringe plunger has been firmly set onto the slider hook and that the syringe has been correctly set and that the flange of the syringe has been inserted into the flange holder. The elevation difference between this product and the patient should be kept to a minimum. Drug solution may not be delivered correctly due to rapid infusion caused by displacement of the syringe plunger from the slider hook, incorrect detection of the syringe size, etc.

Cautions

- When setting a plunger into the slider, check that there is no gap between the plunger of the syringe and the slider. If the plunger is set with gap between them, air or drug solution may be drawn in by the slider hook.
- If the LCD displays the syringe displacement icon and the buzzer sounds, check the syringe has been set correctly, and re-set if necessary. If the syringe is not installed at the correct position, the solution delivery cannot be started.
- Even after the syringe is set in the proper position, if any fault is observed, immediately stop using the product and contact TERUMO trained service technicians.
- When installing a syringe, make sure that the syringe size and brand displayed on the LCD correspond to the syringe. If they do not correspond with each other, it may cause unintended solution delivery.
- When installing a syringe, pull forward the clamp and turn. Then, after the syringe is set, turn the clamp back to the original position and gently push it to fix the syringe. In addition, the slider should be moved while holding the clutch. Any unreasonable operation or excess operation will cause failure of this product.
- When using a syringe (specified sterile syringe for general use) for this product, ensure that the nominal capacity of the syringe is not exceeded. If exceeded, the accuracy of flow rate or the alarm function cannot be guaranteed.

Setting Flow Rate


Set the flow rate per hour.

- Use the Dial to adjust the flow rate in the flow rate screen.

Note

- At this time, the flow rate screen appears in dark blue (selected).

Notes

- The flow rate can be set within the range of the following values.

Setting range:

0.01 to 150.00 mL/h

Note that the upper limit of the flow rate can be changed in the following ranges:

0.01 to 150.00 mL/h (when using syringe of 5 mL)

0.01 to 300.00 mL/h (when using syringe of 10 mL, 20 mL, 30 mL)

0.01 to 1200.00 mL/h (when using syringe of 50/60 mL)

Step:

0.01 mL/h step (0.01 to 10.00 mL/h)

0.10 mL/h step (10.00 to 100.00 mL/h)

1.00 mL/h step (100.00 to 1200.00 mL/h)

- When 150.00 mL/h is reached, the confirmation message "This is the hard upper limit" is displayed on the LCD.

Notes


- For safety, the Dial is designed to make the initial half turn without changing the value on the display. When the value starts to change, the buzzer sounds momentarily for confirmation. When you stop rotating the Dial and start it again in less than 1.5 seconds, the value changes without the initial half turn.
- Rotating the Dial while pressing the Stop switch allows changing the setting at the rate of 100 times.
- When the Dial is turned, a line appears under the numerals to indicate the currently changing digit. (When the Dial is rotated while pressing the Stop switch, an underline appears to the second digit to the left from the normal digit.)
- The set flow rate can be changed during solution delivery. When the Dial is turned on the flow rate screen during solution delivery, the flow rate setting area turns dark blue where the set flow rate increases or decreases in accordance with the Dial operation. The Start switch confirms the set flow rate and changes the solution flow rate. (When the Start switch is not pressed within approx. 5 seconds, it returns to the set flow rate before change.)

Setting VTBI

Set the VTBI dose for the infusion (Default: disabled).

To activate this feature, contact TERUMO trained service technicians.



Setting VTBI to enable (ON) also sets VTBI time to enabled (ON).

1) Press the Display select switch to select the VTBI display.

Notes


- The VTBI tab is selected.
- The flow rate display shows the VTBI display and the VTBI time display.
- The VTBI display appears in dark blue (selected).

2) Use the Dial to adjust the VTBI.

Note

- When the flow rate has already been entered, the VTBI time display indicates the remaining time and VTBI time calculated from the flow rate and VTBI.

Notes


- The VTBI can be set from 0.10 to 9999.00 and "----" (Unlimited).
- When the dose rate is set in mass units, setting in mass units is possible. (VTBI can be set from 0.01 ng to 9999.99 kg and "----" (Unlimited).) (TE-SS800 only)
- If the screen is left non-operated for 10 seconds, it returns to the flow rate screen.
- VTBI can be changed only when delivery is stopped.
- To set to unlimited, rotate half in the direction of decreasing the Dial when VTBI is "0".
- Rotating the Dial while pressing the Stop switch allows changing the setting at the rate of 100 times.

Caution

- In the case of a solution delivery with no setting for the VTBI ("----"; Unlimited), the pump will not stop until the Syringe Empty is issued due to running out drug solution, and therefore the product should be stopped before running out of the drug solution. It is recommended to set a VTBI slightly less than the volume of the drug solution container.

Setting VTBI Time


Set up the time period to complete the dose (Default: disabled).
The flow rate is recalculated from the VTBI and VTBI time specified.
To activate this feature, contact TERUMO trained service technicians.
Setting VTBI to enable (ON) also sets VTBI time to enabled (ON).

- 1) Press the Display select switch to select the VTBI time display.

Note


- The VTBI time display appears in dark blue (selected).

- 2) Use the Dial to adjust the VTBI time.

Notes

- The VTBI time can be set between 00 h 01 min and 99 h 59 min.
- When the VTBI is set to "----" (Unlimited), the VTBI time cannot be set.
- If the screen is left non-operated for 10 seconds, it returns to the flow rate screen.

Cautions

- The flow rate is updated according to the VTBI and VTBI time. Before starting solution delivery, check the value.
- When using the VTBI time setting function, the setting cannot take a flow rate calculation result outside the flow rate setting range, and therefore it should be reconfigured. [----] or [0.00] is displayed, not allowing solution delivery.

◆ Keep Vein Open (KVO) Function

When the volume delivered reaches the VTBI, the KVO function is activated for the prevention of thrombus formation, and thus the flow rate is automatically changed to 0.1 mL/h. If the set flow rate is lower than the flow rate of the KVO, solution delivery continues at the set flow rate. If you wish to change the flow rate of the KVO, contact TERUMO trained service technicians.

Notes

- The KVO is a flow rate set in advance to operate under certain conditions in order to keep the patient line open.
- According to IEC 60601-2-24, clause 2.116, KVO (Keep Vein Open) is commonly used as a synonym for KOR (Keep Open Rate).


◆ Stopping Keep Vein Open

- 1) Press the Back/Mute switch twice to stop the alarm and return to flow rate screen.
- 2) Press the Stop switch to cancel the Keep Vein Open (KVO) function and the operation stops.

Operation Procedure

Priming



Before connecting the infusion line to the patient, perform priming to extract air from the line.

- **Press and hold the Purge switch to perform priming.**

Notes

- During the priming operation, the buzzer sounds and the operation indicator flashes green.
- During the priming operation, the priming volume is displayed on the LCD.
- If any other switch is pressed while the Purge switch is pressed, the safety function detects it as an abnormal operation and stops purging.
- After the priming operation is finished, the screen transits to the flow rate screen.


Notes

- The syringe size and the volume purged are as follows:
 - When using syringe of 5 mL: approx. 150 mL/h
 - When using syringe of 10 mL: approx. 300 mL/h
 - When using syringe of 20 mL: approx. 400 mL/h
 - When using syringe of 30 mL: approx. 500 mL/h
 - When using syringe of 50/60 mL: approx. 1200mL/h
- With the Purge switch pressed and held, the volume delivered display shows the volume purged.
- The priming function will not work when the syringe pump is stopped by an alarm.

Cautions

- Before inserting an intravenous needle or connecting an infusion line with the connecting site, always prime the infusion line by pressing the Purge switch. This product may cause harm to a patient due to air injection, since it does not have a function to detect air in the infusion line.
- Failure by any mechanical gap may cause considerable delay in the start of the infusion. This is especially the case when a large syringe (e.g. 50 ml) is used at very low flow rate settings and/or when anti-siphon valves are used in the infusion line. Use the priming function to remove any mechanical gap: between the syringe plunger and the slider; between the flange and the syringe pump slit.
- When purging the solution, note that the volume purged is not added to the volume delivered. If this is not considered, the volume delivered will differ from the actual infusion volume.
- When purging a drug solution with high viscosity through a thin intravenous needle, an Occlusion alarm is likely to be issued even if the infusion line is not occluded. To deliver a solution with high viscosity, deliver it at the rate of 100 mL/h or lower, not by purging it. If purging is continued, an Occlusion alarm may frequently be issued or the solution may be unable to be delivered.

Inserting a Needle

1) Check


- The pump is stopped.

2) Insert the needle according to the internal hospital procedures.

Warning

- This product does not have an alarm function to detect extravascular infusion due to intravenous needle coming out of vein. Periodically check the puncture site.


Starting Solution Delivery

- Check the flow rate setting and press the Start switch to start solution delivery.

Notes


- The buzzer sounds.
- The operation indicator flashes green.
- The volume delivered increases.

Volume delivered display


Notes

- When the pump is fully set but not started, after 2 minutes the LCD displays the Start Reminder screen (yellow) and the buzzer sounds.
- Use the standby function (see page 40) if you want to temporarily extend the Start Reminder.

When Infusion Cannot Start

The alarm sounds and the following warnings display. (For details, see page 104 "Troubleshooting")


No Flow Rate alarm (yellow)

The flow rate has not been set. Set the flow rate and press the Start switch again.

Cautions

- If soft and hard limits are not activated before starting solution delivery, check the settings for any mistakes such as an error in the digit of flow rate. This product does not have a function to determine whether or not the set values are appropriate, possibly resulting in excessive infusion or insufficient infusion to a patient.
- If any fault is observed, immediately stop the solution delivery and check the setting and the attachment status of the syringe.


Stopping Solution Delivery


1) Press the Stop switch to stop the solution delivery.

Checks


- The operation indicator turns off and the operation stops.
- The volume delivered stops increasing.

2) Close the infusion line.

3) Pull the clamp and then turn.


4) Open the clutch and remove the syringe.

Caution

- When removing the syringe installed in this product, first close the three way stopcock etc. of the infusion line. Otherwise excessive infusion (siphoning free flow by gravity) of the drug solution may occur.


Turning the Power Off

- When the operation ends, press and hold the Power switch for 3 seconds or more to turn the power off.

Notes


- While holding the Power switch, the power off buzzer sounds, then the LCD displays count-down arrows and turns off.

- If the Power switch is pressed and held during solution delivery, a buzzer sounds indicating that solution delivery is in progress. The buzzer continues to sound for approx. 3 seconds until the power is turned off.

Clearing the Volume Delivered

The volume delivered can be cleared.

- 1) Press the Display select switch repeatedly until the volume delivered display is selected.

Note

- The volume delivered display appears in dark blue (selected).



Note

- The volume delivered can be cleared when solution delivery is in progress or when stopped.

- 2) Press the Confirmation switch.

Note

- The buzzer sounds, the volume delivered is cleared to "0", and the screen returns to the flow rate screen.

Notes


- When resuming the solution delivery after the value of the volume delivered is cleared to "0" mL, check the flow rate and the VTBI, and re-set as required.
- To check the volume delivered after clearing the value, use the history function (See page 69).
- Turning the power off clears the volume delivered to "0".
- If the screen is left non-operated for 10 seconds, it returns to the flow rate setting display.

Standby Function

The standby function can extend the Start Reminder.

Note

- The syringe is set and the solution delivery is ready to be started.

- Press and hold (2 seconds or more) the Back/Mute switch. (The standby mode is entered.)**

Notes


- The operation indicator flashes green and red alternately.
- The LCD displays the message "Stand By".

Notes

- The default standby time is 24 hours. To change this, contact TERUMO trained service technicians (the time can be changed to 24 hours or less).
- Turning the Dial or pressing any switch cancels standby mode. (Pressing the Start or Purge switch starts delivery of solution or starts the Purge function.)
- The standby mode cannot be entered when an alarm has been issued. Eliminate the cause of the alarm first.

Replacing Syringe to Continue Solution Delivery

- 1) Press the Stop switch to stop the solution delivery.

Note

- For details on stopping solution delivery, see page 37.

- 2) Close the infusion line.

- 3) Remove the syringe.

Note

- For details on removing the syringe, see page 37.

- 4) Check that the syringe brand to be used matches the syringe brand displayed on the LCD, and then place the syringe.

Note

- For details on setting the syringe, see page 27.

- 5) Perform priming before starting solution delivery to the patient, and keep the plunger pushed by the slider.


Note

- For details on priming, see page 34.

- 6) Clear the volume delivered as required.

Note

- For details on clearing the volume delivered, see page 39.


- 7) Switch the solution delivery to the patient side, open the infusion line, and then press the Start switch to start solution delivery.

Checks

- The flow rate is correctly set.
- The operation indicator flashes green.
- The volume delivered increases.

Selecting Dose Mode (Weight Mode) (TE-SS800 only)


The dose unit setting can be switched.

- 1) Press and hold (2 seconds or more) the Display select switch when solution delivery stopped.

Note


- The setting change display appears and then the dose mode select (Dose) display appears.

- 2) Press the Confirmation switch.

Note




- The existing dose mode (light blue) changes to the selected status (dark blue).


- 3) Turn the Dial or press the Display select switch to select "µg/kg/min" or "mg/kg/h".

Notes

- Turning the Dial or pressing the Display select switch moves the white frame.
- The dose mode goes back to "mL/h" when "mL/h" is selected.








- 4) Press the Confirmation switch to confirm the change of dose mode.

Note

- The display for the unit selected changed dark blue.

- 5) When "Change Dose Mode?" is displayed, press the Confirmation switch.
- 6) The dilution and the weight are displayed in order for them to be changed if necessary. For each of these items, turn the Dial to adjust the value and press the Confirmation switch to set.
- 7) When "OK?" is displayed, press the Confirmation switch to confirm all the settings.

Parameter input display


Dilution (Drug Conc.) input screen

Body weight input screen


Parameter input display

Notes

- The screen transition is performed by pressing the Confirmation switch as described. (The parameter input screens are displayed in the order above.)
- When the Back/Mute switch is pressed in the "parameter input display", the screen returns to the previous menu screen.
- When the Back/Mute switch is pressed before inputting the parameters, the screen returns to the "parameter input display".
- When the Back/Mute switch is pressed after inputting the parameters, "Back to the previous page ?" is displayed on the LCD. By pressing the Confirmation switch, parameter input is cancelled, after which the screen returns to the "parameter input display".
- After the setting is completed, press the Display select switch until info tab **i** is selected. The setting information appears (other than mL/h dose mode).
- Dilution and drug concentration have the same meaning.


Other Operation Procedures

8) Turn the Dial to adjust the dosage.


For example: $\mu\text{g/kg/min}$

For example: mg/kg/h

9) Set the VTBI. (Only if the VTBI function is enabled. See page 32 "Setting VTBI")

10) Press the Start switch to start solution delivery.

Notes

- The flow rate is correctly set.
- The operation indicator flashes green.
- The volume delivered increases.

Notes

Calculation Formula of Flow Rate

$$\text{Flow rate (mL/h)} = \frac{\text{Dose rate } (\mu\text{g/kg/min.}) \times \text{Body weight (kg)}}{\text{Dilution (mg/mL)} \times 1000^{*1}} \times 60^{*2}$$

- *1: Coefficient to convert "minute (min.)" to "hour (h)"
- *2: Coefficient to convert "mg" to "μg"

- Based on dilution, the weight and the dosage, the flow rate (mL/h) is automatically calculated, and the flow rate screen shows the result. The dose unit can be set in the range below.

$\mu\text{g/kg/min}$:

Dilution: 0.01 to 999.00 mg/mL
Weight: 0.1 to 300.0 kg
Dosage: 0.01 to 999.00 $\mu\text{g/kg/min}$

mg/kg/h :

Dilution: 0.01 to 999.00 mg/mL
Weight: 0.1 to 300.0 kg
Dosage: 0.01 to 999.00 mg/kg/h

- If the calculation of flow rate is outside of the setting range, the setting is blocked and should be reconfigured.
- Switching the dose unit changes the values of the dilution, dosage and the weight to "0".
- Selecting the same dose unit also changes the values of the dilution, dosage and the weight to "0".
- It is possible to add/remove some of the dose units. (e.g.: $\mu\text{g/kg/min}$ and mg/kg/h can be made unavailable. For details, contact TERUMO trained service technicians.)
- Even after the setting is completed, the unit returns to mL/h when the power is turned off.
- Patient name, Patient ID, Height, Gender and Age can be acquired via communication with an external device. For the setting method, contact TERUMO trained service technicians.

Selecting Dose Mode (Drug Library) (TE-SS800 only)

The dose mode can be selected from the drug library.

Use the dedicated editing software (TERUFUSION Drug Library Manager or TERUFUSION Software Package) for registration. For details, contact TERUMO trained service technicians.

- 1) Check the selected profile. (See page 83 “Profile (TE-SS800 only)” for how to select a profile.)



Note

- Profiles allow drugs to be classified according to their use and registered in the pump. Selecting a profile can limit the available libraries.

- 2) Press and hold (2 seconds or more) the Display select switch when solution delivery stopped.

Note


- The setting change display appears and then the dose mode select (Dose) display appears.

- 3) Press the Confirmation switch.

Note



- The existing dose mode (light blue) changes to the selected status (dark blue).


- 4) Turn the Dial or press the Display select switch to select the Library.

Notes

- The item, “Library”, changes has changed from (light blue) to the selected status (dark blue).
- Turning the Dial or pressing the Display select switch moves the white frame.



5) Press the Confirmation switch.

Notes

- The screen changes to the library selection display.
- The item, "Library", changes from (light blue) to the selected status (dark blue).

Screen example: When "Profile1" is selected, "Library1", "Library2", "Library3", "Library4" or "Library5" can be selected.




Note


- The selected profile name is displayed between "Library" selected and the drug name in the dose mode.
- (See page 83 "Profile (TE-SS800 only)" for how to select a profile.)

6) Turn the Dial or press the Display select switch to select the library.

Notes

- The previous setting has been selected (dark blue).
- Turning the Dial or pressing the Display select switch moves the white frame.



7) Press the Confirmation switch.

Notes

- The screen changes to display the drug information of the selected library such as the dosage rate setting range and dose unit.

- The detailed information can be displayed by turning the Dial or pressing the Display select switch.


- See the dedicated editing software (TERUFUSION Drug Library Manager or TERUFUSION Software Package) for how to set detailed information such as the dosage rate setting range and dose unit.

8) Check the detailed information and press the Confirmation switch.


Notes

- Clinical advisory is displayed. (* if clinical advisory is set)


- Clinical advisory is a function to display a cautionary message about the selected drug.

See the dedicated editing software (TERUFUSION Drug Library Manager or TERUFUSION Software Package) for how to set clinical advisory for each drug.



Other Operation Procedures


- 9) Press the Confirmation switch until “Change Dose Mode?” is displayed on the LCD and press the Confirmation switch again.

Note

- The display for the Library selected appears.


- 10) The dilution and the weight are displayed in order to be changed if necessary. For each of these items, turn the Dial to adjust the value and press the Confirmation switch to set.

Notes


- The screen transition is to be performed by pressing the Confirmation switch as described. (the parameter input screens are to be displayed in sequence)
- When the Back/Mute switch is pressed in the “parameter input display”, the screen returns to the previous menu screen.
- When the Back/Mute switch is pressed before inputting the parameters, the screen returns to the “parameter input display”. When the Back/Mute switch is pressed after inputting the parameters, “Back to the previous page ?” is displayed on the LCD. By pressing the Confirmation switch, parameter input is cancelled, after which the screen returns to the “parameter input display”.
- Setting dilution is possible only when the default dilution of the library is 0.00 (not set).
- The concentration value cannot be manually changed if the default value is specified in the library.
- The ranges of the dilution and weight are as follows:
Dilution: 0.01 to 999.00
Weight: 0.1 to 300.0 kg
- The dilution and weight specified are reset when the dose unit/library has been changed, or the power has been turned off.


- 11) Select “OK?” and press the Confirmation switch.

Note

- The flow rate setting display is shown.

Library display



- Colour tag
Selection is possible from 16 colours.
See the dedicated editing software (TERUFUSION Drug Library Manager or TERUFUSION Software Package) for how to set an advanced dose mode and colour tag.
- Advanced dose mode icon
The icon to indicate the advanced dose mode is displayed.

 Advanced dose mode icon	The icon is displayed when the advanced dose modes below can be selected. For details of each dose method, see "TERUFUSION Drug Library Manager" or "TERUFUSION Software Package". Note • If no advanced dose mode can be set for the selected library, the advanced dose mode icon is not displayed.	
	 Interval/Intermittent mode icon	Two specified dosage rates are alternately delivered at regular intervals for the specified number of steps.
	 Multi-step mode icon	The speed of solution delivery increases or decreases at the specified levels.
	 Delayed Start mode icon	Solution delivery starts after the specified time has elapsed.
 Combination dose mode icon		For details, contact TERUMO trained service technicians.

Notes

- After the programming is completed, press the Display select switch until info tab  is selected.
The programming information appears.
- The "Combination Dose Mode" which can be used in combination with a rack system (optional accessory) is available in the library. Contact TERUMO trained service technicians for details.
- "Pressure alarm" is available in the library as a safety function.
See "TERUFUSION Drug Library Manager or TERUFUSION Software Package" or contact TERUMO trained service technicians for details.
- "Come and See ME" is available in the library as a special function. See "TERUFUSION Drug Library Manager or TERUFUSION Software Package" or contact TERUMO trained service technicians for details.


12) Turn the Dial to adjust the dosage.

Notes

- The calculation result of flow rate is within the flow rate setting range.

Dosage rate

- The range of the dosage rate is as follows:
Dosage: 0.01 to 999.00
- The flow rate (mL/h) is automatically calculated from the specified dilution, weight and dosage rate, and displayed in the flow rate display.
- If the flow rate reaches the soft limit of the lower limit or the soft limit of the upper limit, the "Soft lower limit reached, Override?" or "Soft upper limit reached, Override?" confirmation message is displayed on the LCD.
- If the flow rate reaches the hard limit of the lower limit or the hard limit of the upper limit, the "This is the hard lower limit" or "This is the hard upper limit" confirmation message is displayed on the LCD.

Caution


- When using the drug library mode (TE-SS800 only), the setting cannot take a flow rate calculation result outside the flow rate setting range, and therefore it should be reconfigured. [----.--] or [0.00] is displayed, not allowing solution delivery.

13) Set the VTBI. (Only if the VTBI function is enabled. See page 32 "Setting VTBI")

14) Press the Start switch to start solution delivery.

Notes

- The flow rate is correctly set.
- The operation indicator flashes green.
- The volume delivered increases.
- Even after the setting is completed, the unit returns to mL/h when the power is turned off.

Setting Special Functions

The following settings are enabled by changing the internal settings of this product. If necessary, contact TERUMO trained service technicians.

Setting item	Description of function
VTBI setting	Allows you to set the VTBI. (For details, see page 32 "Setting VTBI".)
VTBI time setting	Allows you to set the time limit. (For details, see page 33 "Setting VTBI Time".)
Setting step 100 times function	When turning the Dial while pressing the Stop switch, the plus minus rate of the set value (flow rate etc) accelerates. (For details, see page 5 "Operation Panel".)
Syringe brand setting	The syringe brand can be selected after the power is turned on and the syringe installation has finished.
Flow rate range setting function A (Soft limits) Flow rate range setting function B (Hard limits)	Sets the limitation of the flow rate in 2 levels; the soft limit (can be exceeded the upper limit after confirmation) and the hard limit (cannot be exceeded). • The setting step 100 times function only after the flow rate exceeds the soft limit. • When setting the flow rate over the soft limit, the flow rate will be displayed in orange. • When the drug library mode or body weight dose mode ($\mu\text{g/kg/min}$, mg/kg/h) is selected, the flow rate range setting function is not available. (See page 42 "Selecting Dose Mode (Weight Mode) (TE-SS800 only)".)
Change the Nearly Empty alarm setting	Allows you to set the Nearly Empty alarm position based on the time before the syringe is empty.
Change the KVO flow rate (Only if the VTBI function is enabled)	Allows you to set the KVO flow rate within the range of 0.10 to 3.00mL/h.
Hands Free Bolus setting	When bolus volume and bolus time are set, flow rate is calculated and the bolus infusion is available. (For details, see page 61 "Hands Free Bolus".)
Volume delivered unit select function (TE-SS800 only)	In the case of mass system dose mode, the unit for volume delivered or VTBI can be changed from volume (mL) to mass (g). (TE-SS800 only)
Alarm melody select function	Select between standard specified melody and TERUMO original melody.
Night mode setting	If the device is left non-operated for 30 seconds between 21:00 to 5:00 (customizable), the brightness can be automatically lowered by two levels (customizable).
Change the standby duration	Up to 24 hours with 10 minute intervals.
Change the startup dose mode (TE-SS800 only)	The dose mode at power ON can be selected. mL/h, $\mu\text{g/kg/min}$, mg/kg/h and Library can be selected. (TE-SS800 only)
" μ " display switching function (TE-SS800 only)	Switches the dose unit from " μ " to "mc".

Changing Settings on the Menu

See pages 10 and 11 for the menu list.

Notes

- If you press the Start switch whilst changing the settings from the menu, solution delivery starts (when the syringe has been attached, the flow rate, etc., have been set, and the product is ready for operation).
- Occlusion detection pressure, keypad lock, bolus, advanced dose mode and combination dose mode on Menu 1 can be selected or set during solution delivery. Other items are not available.
- How to exit menus
 - Method 1) Press and hold the Display select switch in the menu for 2 seconds or more.
If the setting has not been completed, after the confirmation screen is displayed, cancel the setting to go to the flow rate screen.
 - Method 2) Press the Back/Mute switch repeatedly.
However, in the setting select screen (deepest layer), select a setting and go back to the previous screen.
- If an alarm occurs before pressing the Confirmation switch or the Back/Mute switch while setting the menu, any changes will return to the settings before the alarm was issued.
- Time before the screen returns to the flow rate screen.
 - (1) Level in which menu tab has been selected:
Goes back to the flow rate screen after 10 seconds of no operation.
 - (2) Levels on and after which menu tab has been specified:
Goes back to the flow rate screen after 2 minutes of no operation.

Occlusion Detection Pressure

The occlusion detection pressure can be set in 10 levels.


Setting	Occlusion detection pressure	Setting	Occlusion detection pressure
1	± 10 kPa	6	± 70 kPa (Default)
2	± 20 kPa	7	± 90 kPa
3	± 30 kPa	8	± 100 kPa
4	± 40 kPa	9	± 110 kPa
5	± 50 kPa	10	± 120 kPa

Warning

- After the occlusion detection pressure is selected, check that the setting has been adjusted correctly.

Notes

- If using on "Level 10", the pressure of an infusion line increases at occlusion, possibly resulting in the syringe connection etc. being disconnected. Periodically check that the connection of the infusion line is secured.
- If using on "Level 1", an Occlusion alarm is frequently issued, especially when using a drug solution with high viscosity, possibly resulting in infusion being interrupted. In such cases, review the set value of the occlusion pressure by using the graph of the history of occlusion pressure value as a guide.
- When an Occlusion alarm is issued, the occlusion icon  may not reach the orange level (more than 75%) for the following reasons:
 - When the slider pull-back function is activated in order to mitigate the occlusion status
 - When the syringe is fully pressed (stopped after the Nearly Empty alarm or an internal pressure equivalent to the set "Level 4" occlusion pressure value is detected)
- After an Occlusion alarm is issued internal pressure may build up even if the occlusion icon has not turned orange.
- Use the occlusion icon  as follows:
 - To check the in-line pressure when using a small syringe and a small diameter administration set and/or infusing a highly viscous solution.
 - To make sure the infusion line has been opened.


- Press and hold the Display select switch for 2 seconds or more to display the menu 1 screen.

Note

- The menu 1 screen is displayed and the icons are ready to be selected.


Changing Settings on the Menu

- 2) Turn the Dial or press the Display select switch to select the occlusion detection pressure tab .

Note

- The occlusion detection pressure setting change display has appeared.




Note

- In the graph, the horizontal axis indicates the time and the vertical axis indicates the occlusion detection pressure.

- 3) Press the Confirmation switch.

Note


- The existing occlusion detection pressure (light blue) changes to the selected status (dark blue).

- 4) Turn the Dial or press the Display select switch to select occlusion detection pressure.


Note

- Turning the Dial or pressing the Display select switch moves the white frame.

Notes

- Upper half: Occlusion detection pressure setting selection
- Lower half: Occlusion history graph
 - In the graph, the horizontal axis indicates the time and the vertical axis indicates the occlusion detection pressure.
 - Displays a graph of the history of occlusion pressure value in the past 2 hours. The horizontal line in the graph indicates the currently selected occlusion detection pressure. (This history is deleted when the power is turned off.) The graph is updated every 2 minutes.

5) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

-  : Confirm and go back to the flow rate screen.
-  : Confirm and go back to the previous screen.

Notes

- To cancel editing, select the original value and press the Back/Mute switch to return to the menu 1 screen.
- The occlusion icon  in the status icon display shows the selected occlusion detection pressure.


Changing Settings on the Menu

Keypad Lock Function

The panel operation can be disabled to prevent inadvertent operation while solution delivery is stopped or in-progress.

Notes



- Even if the keypad lock is enabled, the Power switch and the Start switch can be operated. The keypad lock release operation (press and hold the Back/Mute switch) is also allowed when keypad lock is enabled.
- Even if the Stop switch is pressed, the solution delivery cannot be stopped when keypad lock is enabled. To stop solution delivery, first release keypad lock.

- 1) Press and hold the Display select switch for 2 seconds or more to display the menu 1 screen.

Note



- The menu 1 screen is displayed and the icons are ready to be selected.

- 2) Turn the Dial or press the Display select switch to select the keypad lock tab.

Note


- The keypad lock ON/OFF setting display has appeared.


- 3) Press the Confirmation switch.

Note

- The existing setting (light blue) changes to the selected status (dark blue).


4) Turn the Dial or press the Display select switch to select ON.

Note

- Turning the Dial or pressing the Display select switch moves the white frame.

5) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

: Confirm and go back to the flow rate screen.

: Confirm and go back to the previous screen.

Note


- To cancel editing, select the original value and press the Back/Mute switch to return to the menu 1 screen.

Checks

- The keypad lock icon  appears in the status icon display.
- The keypad lock message appears when the Dial or a switch (other than the Power switch and Start switch) is operated.

◆ To release keypad lock

- Press and hold the Back/Mute switch for 2 seconds or more to release the keypad lock.

Check

- The keypad lock icon has disappeared.

Note

- Even while the keypad lock is on, the Power switch and the Start switch are available.


Note

- In the case of a high-priority alarm, the keypad lock function is automatically switched off. However, keypad lock will not be disabled by Start Reminder.
(If the status allows startup, the solution delivery can be started while keypad lock is enabled.)

Bolus

◆ Hands On Bolus

There are two ways for selecting bolus either through the menu or pressing the Purge switch while pressing the Back/Mute switch. In addition, every time the Purge switch is pressed while Back/Mute switch is pressed, you can switch between the Hands On Bolus and Hands Free Bolus.



- 1) Press and hold the Display select switch for 2 seconds or more to display the menu 1 screen.

Note

- The menu 1 screen is displayed and the icons are ready to be selected.



- 2) Turn the Dial or press the Display select switch to select the bolus tab.


- 3) Press the Confirmation switch.

Note

- The existing setting (light blue) changes to the selected status (dark blue).


Changing Settings on the Menu

4) Press the Confirmation switch to confirm the selection.

Note


- The current bolus flow rate setting screen appears.

5) Turn the Dial to set the bolus flow rate.

Note




- The bolus flow rate can be set by 100 mL/h in the following value range.
 - When using syringe of 5 mL: 100 mL/h
 - When using syringe of 10, 20, 30 mL: 100 mL/h to 300 mL/h
 - When using syringe of 50/60 mL: 100 mL/h to 1200 mL/h


6) Press and hold the Purge switch to start the bolus.

Notes

- The operation indicator flashes green.
- The bolus volume delivered increases.
- The alarm melody sounds during solution delivery.

Solution volume delivered

7) Release the Purge switch to stop the bolus.


Notes

- The screen returns to the flow rate screen after the bolus stops.
- If the Purge switch is pressed again within 1.5 seconds after being released, Hands On Bolus continues and the bolus volume delivered display is maintained. If exceeding 1.5 seconds after the Purge switch has been released, the screen returns to the flow rate screen.
- While the Purge switch is pressed, the volume delivered display shows the bolus volume delivered. When releasing the Purge switch, the display shows the sum of the solution volume delivered. (The bolus volume delivered is added to the sum of the solution volume delivered.)
- When bolus infusion during solution delivery is completed, solution delivery starts at the same flow rate as before the bolus infusion.
- The dose setting (flow rate) of Hands On Bolus remains set, and is used when the next Hands On Bolus is performed, until the power is turned off or the dose mode is switched.

◆ Hands Free Bolus

Setting the Bolus VTBI and Bolus VTBI time enables automatic calculation of the Bolus flow rate (mL/h) and execution of bolus infusion without holding down the Purge switch. (Hands Free Bolus function)

To set the Hands Free Bolus function to ON, contact TERUMO trained service technicians.



1) Press and hold the Display select switch for 2 seconds or more to display the menu 1 screen.

Notes


- The menu 1 screen is displayed and the icons are ready to be selected.
- VTBI needs to be set before setting this function if VTBI has been enabled.

2) Turn the Dial or press the Display select switch to select the bolus tab.



Changing Settings on the Menu



3) Press the Confirmation switch.

Note

- The existing setting (light blue) changes to the selected status (dark blue).


4) Turn the Dial or press the Display select switch to select the Hands Free Bolus.

5) Press the Confirmation switch.



Notes

- The bolus setting display appears.


- When the flow rate screen is displayed, the Hands Free Bolus screen can be displayed by pressing and holding the Back/Mute switch and pressing the Purge switch twice.

6) Turn the Dial to set the bolus VTBI.




Note

- The Bolus VTBI can be set within the range below.
VTBI: 0.01 to 999.00


- 7) Press the Display select switch to select the Bolus VTBI time or the Bolus flow rate.

Bolus VTBI time

Bolus flow rate

- 8) Turn the Dial to set the Bolus VTBI time or the Bolus flow rate.

Notes

- The Bolus VTBI time or the Bolus flow rate (mL/h) has been calculated from the bolus VTBI, and is displayed.
- The start up icon is displayed.

Ready for bolus icon

Bolus VTBI time


Bolus flow rate

- The Bolus flow rate can be set by 100 mL/h in the following value range.
 - When using syringe of 5 mL: 100 mL/h
 - When using syringe of 10, 20, 30 mL: 100 mL/h to 300 mL/h
 - When using syringe of 50/60 mL: 100 mL/h to 1200 mL/h
- The Bolus VTBI time can be set within the range below.
Bolus VTBI time: 00min01s to 60min00s

Caution

- When using the Hands Free Bolus function, the setting cannot take a flow rate calculation result outside the flow rate setting range, and therefore it should be reconfigured.
[----.--] or [0.00] is displayed, not allowing solution delivery.

Changing Settings on the Menu

9) Press the Purge switch to start the bolus.

Notes

- The operation indicator flashes green.
- The bolus volume delivered increases.
- The alarm melody sounds during solution delivery.

10) The bolus ends when the VTBI is reached.

Check




- The screen returns to the flow rate screen when the VTBI has been reached and the bolus has been completed.

Notes

- The volume delivered display area displays the amount of bolus during bolus infusion, and displays the sum of the amount of bolus and the volume delivered when the VTBI is attained and the bolus ends.
- When bolus infusion during solution delivery is completed, solution delivery continues at the flow rate before the bolus infusion.
- Pressing the Stop switch forcibly stops Hands Free Bolus.

Advanced Dose Mode (TE-SS800 only)

An advanced dose mode can be selected from the three advanced dose modes in the table below according to the selected drug. (Before selecting the dose mode, it is necessary to register a library for which an advanced dose mode can be selected and to select the library.) See page 45 "Selecting Dose Mode (Drug Library) (TE-SS800 only)" for how to select the drug.

Advanced dose mode	Details
 Interval/Intermittent Mode	Two specified dosage rates are alternately delivered at regular intervals for the specified number of steps.  Input parameters 1) Dosage rate 1*¹ 2) Dosage rate 2 (Displays: Final Flow Rate) 3) 1 cycle time (Displays: Time/Cycle) 4) Time interval (Displays: Time) 5) Number of steps (Displays: Number of Steps) { If Dose Unit = mL/h : 0.00 to 1200.00*² { If Dose Unit = other than mL/h : 0.00 to 999.00*² 00h01min to 23h59min 00h01min to 23h59min 1 to 255 steps
 Multi-step Mode	The speed of solution delivery increases or decreases at the specified levels (up to 10 levels).  Input parameters 1) Dosage rate 1*¹ 2) Dosage rate 2 (Displays: Final Flow Rate) 3) Time (Displays: Time to Final Flow Rate) 4) Number of steps (Displays: Steps) { If Dose Unit = mL/h : 0.00 to 1200.00*² { If Dose Unit = other than mL/h : 0.00 to 999.00*² 00h01min to 23h59min 1 to 10 steps
 Delayed Start Mode	Solution delivery starts after the specified time has elapsed.  Input parameters 1) Dosage rate 1*¹ 2) Dosage rate 2 (Displays: Final Flow Rate) 3) Delay Time (Displays: Delay Time) { If Dose Unit = mL/h : 0.00 to 1200.00*² { If Dose Unit = other than mL/h : 0.00 to 999.00*² 00h01min to 23h59min

*¹ "Dosage rate 1" can be set in the dose rate in the flow rate setting screen.

*² There is no rule stating which of the volumes between dosage rate 1 and dosage rate 2 should be larger.

If both the dosage rates 1 and 2 are not set, solution delivery does not start.

The calculation result of flow rate (mL/h) must be set within the flow rate setting range.

The dosage rate must be set within the setting range for each library.

Note

- When registering Library, determine the range of the dosage rate considering the drug and the dose method. (Depending on the setting of the dosage rate, the occlusion may occur due to the thrombus formation within the infusion line.)

Changing Settings on the Menu


Advanced dose mode icon 1

Dosage rate 1

- 1) Follow step 1) to 11) on "Drug Library" to select the library and set the dosage rate 1.

Notes



- See page 45 "Selecting Dose Mode (Drug Library) (TE-SS800 only)" for how to select the library.
- If the selected library has been enabled to accept an advanced dose, the advanced dose mode icon  is displayed on the LCD. If an advanced dose is unacceptable, no icon is displayed.

- 2) Press and hold the Display select switch for 2 seconds or more to display the menu 1 screen.

- 3) Turn the Dial or press the Display select switch to select the Advanced dose mode tab .




Advanced dose mode tab

- 4) Press the Confirmation switch.

Notes

- The existing setting (light blue) changes to the selected status (dark blue).
- The available advanced dose modes are displayed.



- 5) Turn the Dial or press the Display select switch to select the advanced dose mode.

Note

- Turning the Dial or pressing the Display select switch moves the white frame.


- 6) Press the Confirmation switch.


Check


- The screen changes to the parameter input display.

Interval/Intermittent Mode


Multi-step Mode


Delayed Start Mode

- 7) Turn the Dial to input the parameters.


Changing Settings on the Menu

E. g.) Interval/Intermittent Mode

Notes



- The screen transition is performed by pressing the Confirmation switch as described.
(The parameter input screens are displayed in the order above.)
- When the Back/Mute switch is pressed in the "parameter input display", the screen returns to the previous menu screen.
- When the Back/Mute switch is pressed before inputting the parameters, the screen returns to the "parameter input display".
When the Back/Mute switch is pressed after inputting the parameters, "Back to the previous page ?" is displayed on the LCD.
By pressing the Confirmation switch, parameter input is cancelled, after which the screen returns to the "parameter input display".

8) Select "OK?" and press the Confirmation switch.

Checks


- The input parameters are correct.
- The screen returns to the flow rate screen, and the icon of the selected dose mode is displayed in the lower left of the flow rate display.

9) Press the Start switch to start solution delivery.

Notes

- The solution delivery has started and is executed at the specified dose settings.
- The operation indicator flashes green.
- The volume delivered increases.


History Function

With the history function, operation history (such as date, time, set flow rate, when the solution delivery was started or stopped, or when an alarm occurred) can be checked. The latest records of max. 10,000 items are stored. If 10,000 items are exceeded, the records are deleted from the oldest first.

Notes



- The battery charge level display indicates the current status regardless of the displayed operation history.
- If the pump is dropped then a shock record will be created in the history. This record will be shown at the top of history.

- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note


- The menu 2 screen is displayed and the icons are ready to be selected.

- 2) Turn the Dial or press the Display select switch to select the history tab.

Note


- The history display has appeared.



- 3) Press the Confirmation switch.

Note

- The existing setting (light blue) changes to the selected status (dark blue).

Changing Settings on the Menu

- 4) Turn the Dial or press the Display select switch to select the desired operation history number.

Note

- Turning the Dial or pressing the Display select switch moves the white frame.




- 5) Press the Confirmation switch to display the detailed history and use the Display select switch or Dial to select a history.

Notes

- The operation histories are classified by operation status and identified with the colours at the bottom area of LCD as shown in the table below.
- Screen transition
 - As described in Step 4), the Display select switch and Dial can be used to move up and down in the list.
 - Pressing the Confirmation switch in the list display displays the detailed history.
 - Pressing the Back/Mute switch in the detailed history returns to the list display.
 - As described in Step 5), the Display select switch and Dial can be used to move the previous or next detailed history.

Example of operation history	Colour of bottom area of LCD	Operation status
Info. tab 2012/07/01 01:25 White	White	Operation history information
2012/07/01 09:15 Green	Green	Operation history during solution delivery

Example of operation history	Colour of bottom area of LCD	Operation status
 Pink	Pink	Operation history during stop
 Blue	Blue	Operation history during purge or bolus (Hands On Bolus, Hands Free Bolus)
 Purple	Purple	Operation history during standby

6) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.

Note


- Set the correct time and date by following the procedure described on page 79 prior to use.

Changing Settings on the Menu

Changing the Sound Volume

The alarm sound volume, the operation sound volume, and the Cable Insertion/Extraction sound volume can be set on 4 levels.


Type	Volume setting	Default setting	Note
Alarm Volume	1 to 3, Step	2	0 (very low), 1 (low), 2 (medium), 3 (high), Step (1 → 2 → 3 is repeated)
Key Operating Volume	0 to 3	1	
Cable Insertion Alarm	0 to 3	0	



- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note

- The menu 2 screen is displayed and the icons are ready to be selected.


- 2) Turn the Dial or press the Display select switch to select the sound volume tab.


- 3) Press the Confirmation switch.


Note

- The white frame is displayed.



- 4) Turn the Dial or press the Display select switch to select an item to change.

- 5) Press the Confirmation switch to confirm the selection.

Note


- The display for the item to be changed has appeared.


- 6) Turn the Dial or press the Display select switch to select the volume.

Notes


- Turning the Dial or pressing the Display select switch moves the white frame.
- When the cursor is placed on a desired level, the sound is audible.
- If "Step" for the alarm volume is selected, all 3 levels of volume are demonstrated.

Alarm Volume


Key Operating Volume

Cable Insertion Alarm

Changing Settings on the Menu



- 7) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

: Confirm and go back to the previous screen.

Check

- The selected volume is displayed.



- 8) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.

Changing the Brightness

The brightness can be set in 5 levels. (Default: 5)

By default, the night mode is activated. (For details, see page 51 "Setting Special Functions")



- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note

- The menu 2 screen is displayed and the icons are ready to be selected.



- 2) Turn the Dial or press the Display select switch to select the brightness tab.


- 3) Press the Confirmation switch.

Note

- The existing setting (light blue) changes to the selected status (dark blue).



Changing Settings on the Menu

- 4) Turn the Dial or press the Display select switch to select a brightness level.

Notes

- Turning the Dial or pressing the Display select switch moves the white frame.
- When the cursor is placed on a desired brightness level, the brightness level is shown.



- 5) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

: Confirm and go back to the previous screen.

Checks

- The LCD is displayed with the brightness selected.
- The setting is selected (light blue).


- 6) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.

Maintenance Timer Function

The time period until the maintenance timer display is lit can be set to value between OFF (None) and 36 months in one month increments. The maintenance timer icon is lit when the set period has elapsed. Perform maintenance, and then set the time period again. When the next setting (update) is made after the maintenance icon  is once lit, make the setting again.



Display	Detail
OFF (Default)	No setting
1 to 36	1 to 36 months



- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note

- The menu 2 screen is displayed and the icons are ready to be selected.


- 2) Turn the Dial or press the Display select switch to select the maintenance timer tab .


- 3) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

Note

- The existing setting (light blue) changes to the selected status (dark blue).


Changing Settings on the Menu


- 4) Turn the Dial or press the Display select switch to select the period.

Note

- Turning the Dial or pressing the Display select switch moves the white frame.

- 5) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

: Confirm and go back to the previous screen.

Check

- The date after the set period is displayed.



YYYY/MM/DD

The next maintenance timer alert date is displayed in (YYYY: year, MM: month, DD: date) format.

- 6) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.


Note

- The maintenance icon is lit when the maintenance timer set period has elapsed.

Setting the Date and Time

The date and time can be set. Set YYYY: year, MM: month, DD: date, hh: hour, and mm: minute.



- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note

- The menu 2 screen is displayed and the icons are ready to be selected.

- 2) Turn the Dial or press the Display select switch to select the date and time tab .



Date and time tab


- 3) Press the Confirmation switch.

Note


- The white frame is displayed.



Changing Settings on the Menu


4) Turn the Dial or press the Display select switch to select an item.

5) Press the Confirmation switch.


6) Turn the Dial to change the value.

7) Press the Confirmation switch.


Note
• The display for the next item appears.


8) To continue to change the next item, repeat Steps 6) and 7).

Note

- Press the Confirmation switch to cycle through displays in the order of year, month, date, hour, and then minute.

9) Press the Confirmation switch to change the screen to waiting for a confirmation status (OK?) after the minute setting is changed.

Changing Settings on the Menu


◆ When changing the settings

- 10) Press the Confirmation switch to set a specified date and time with no specified seconds (Default 00 s).


◆ If any modification is required

- 10) Press the Back/Mute switch and when “Back to the previous page ?” is displayed, press the Back/Mute switch again. Return to Step 4) and turn the Dial to select an item to change.


◆ To cancel the settings

- 10) Press the Back/Mute switch and when “Back to the previous page ?” is displayed, press the Confirmation switch to return to the previous screen without confirming the change.

Profile (TE-SS800 only)



Profiles allow drugs to be classified according to their use and registered in the pump. Selecting a profile can limit the drugs available to each individual profile. Up to 30 profiles can be registered. See page 45 “Selecting Dose Mode (Drug Library) (TE-SS800 only)” for more information about drugs.


- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note

- The menu 2 screen is displayed and the icons are ready to be selected.

- 2) Turn the Dial or press the Display select switch to select the profile tab. If no library or profile is registered, the profile tab is not displayed.


- 3) Press the Confirmation switch.

Note

- The existing setting (light blue) changes to the selected status (dark blue).




Changing Settings on the Menu


- 4) Turn the Dial or press the Display select switch to select a profile.

Note

- Turning the Dial or pressing the Display select switch moves the white frame.

- 5) Press the Confirmation switch or the Back/Mute switch to confirm the selection.

: Confirm and go back to the previous screen.


Check

- In "Drug library", only the libraries set in the selected profile are available. See page 45 "Selecting Dose Mode (Drug Library) (TE-SS800 only)" for details on the library.

Note

- See the dedicated editing software (Dose Analyzer with Drug Library or Drug Library Manager) for profile setting and registration of each library.

- 6) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.


Communication Environment (TE-SS800 only)

Wireless LAN can be turned ON/OFF and the communication environment can be set. Up to 10 communication environments can be registered.

◆ Communication environment setting



For details about how to register the communication environment setting to the pump unit, contact TERUMO trained service technicians.



- 1) With the Back/Mute switch pressed and held, press and hold (2 seconds or more) the Display select switch to display the menu 2 screen.

Note

- The menu 2 screen is displayed and the icons are ready to be selected.


- 2) Turn the Dial or press the Display select switch to select the communication environment tab ➡ .


- 3) Press the Confirmation switch.

Note


- The status changes from the icon selection (light blue) to the item selection (white frame).


Changing Settings on the Menu

4) Check that “Communication Mode Select” is selected, and then press the Confirmation switch.


5) Turn the Dial or press the Display select switch to select the communication environment you want to set from the list of registered communication environment.


Notes

- Turning the Dial or pressing the Display select switch moves the white frame.
- The specified communication environments are restricted depending on the set profile.


6) Press the Confirmation switch or Back/Mute switch to confirm the selection.

Check

- The name of the set communication environment is displayed at the upper right of the communication environment setting screen.



7) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.


◆ Wireless LAN ON/OFF setting

- 1) Follow Steps 1) to 3) of “Communication environment setting” (page 85) to display the communication environment setting screen.
- 2) Turn the Dial or press the Display select switch to select Wireless LAN.

Note
• Turning the Dial or pressing the Display select switch moves the white frame.



- 3) Press the Confirmation switch.


- 4) Turn the Dial or press the Display select switch to select ON/OFF.


Notes
• The previous setting has been selected (dark blue).
• Turning the Dial or pressing the Display select switch moves the white frame.

Changing Settings on the Menu



5) Press the Confirmation switch or Back/Mute switch to confirm the selection.

Check

- The selected wireless LAN setting is displayed.

6) Press and hold (2 seconds or more) the Display select switch to return to the flow rate screen.

Overview

Use infrared communication (IrDA), external communication (RS-232C) or wireless LAN to output the status of the pump to a network system (e.g. monitor system or CIS) (The status of the pump cannot be changed from a network system)
Ensure that the network system and the pump use the same communication conditions.

List of communication specifications

Communication method	Infrared communication (IrDA) (TE-SS800 only)	External communication (RS-232C) (TE-SS702 only)	Wireless LAN (TE-SS800 only)
Transmission rate	115200bps	115200bps	1Mbps to 11Mbps
Data length	8 bits	8 bits	-
Start bit	1 bit	1 bit	-
Stop bit	2 bits	1 bit	-
Parity check	None	None	-
Frequency (Europe)	-	-	2412 MHz (ch1) to 2472 MHz (ch13)
Encryption scheme	-	-	WPA2 (AES/TKIP), WPA (AES/TKIP), WEP
Operational mode	-	-	Infrastructure mode

Cautions

- When this product is connected with other devices, particularly note the following:
 - When connecting this product with other medical equipment or network system, check its conformance to IEC 60601-1-1: 2000 and IEC 60601-1-2: 2001 (Amd.1: 2004) prior to use in order to ensure the system safety.
 - When connecting this product with a network system or other equipment, for safety, checking the manufacturer's specifications of the equipment is recommended. In addition, use EMI compliant products for the connection cable.
 - When using the external communication function, pay particular attention, as it is more susceptible to the effect from a radio knife, mobile phone, radio device, defibrillator, etc. Regularly check that this product is operating normally.
 - If connecting with a nurse call, always check in advance with the electrical contractor for compatibility, and connect with a system for which the connection function has been confirmed.
[The functions regarding connection are not guaranteed.]
 - If connecting with a buzzer or lamp instead of a nurse call, it should be used with a capacity of 12V DC, 1A or less.
 - When using the wireless LAN communication function, the effect on other equipment by radio interference should be considered. [It may affect other equipment.]
 - When connecting this product with a network, the settings of this product and the network system need to be compatible. Contact a vendor specialist for the correct settings and allow only the system administrator to perform settings. [Failure to connect with correct settings may result in interfering with the original functionality or performance of this product and affecting the network system.]

<Cautions regarding connecting communication cables, communication devices and the nurse call button>

- When using this product's external communication function to connect a network system and a nurse call button, check the specifications of the monitor manufacturer, the communication cable manufacturer, and the nurse call button manufacturer, and further, check safe operation. Further, when connecting to a nurse call button, do not exceed this product's rating (maximum 1 A). Also, in order to check the safety of the connections as a system, check their compatibility with IEC 60601-1-1: 2000 and IEC 60601-1-2: 2001 (Amd.1: 2004). Not following these cautions could result in interference with other medical electrical equipment or breakdown.

Notes

- Operating time using the internal battery is decreased when using external communication.
- The specified message (free message function) can be displayed on the pump unit via communication. (For details, contact TERUMO trained service technicians.)

Network System Connection

Procedure for Connecting the Network System

Caution

- When connecting this product with a network system or other equipment, for safety, checking the manufacturer's specifications of the equipment is recommended. In addition, use EMI compliant products for the connection cable.


< IrDA (For infrared communication)>

Fix the pump into the optional TERUFUSION Communication Rack System and then turn the power on for both so that the pump and rack system automatically communicate with each other. See the instruction manual for the rack system for communication between rack system and external device.

< RS-232C (For external communication terminal)>

- 1) Check that all the devices are powered off.
- 2) Connect the cable.
- 3) Turn on the devices.

Cable connection example (For RS-232C)


Note

- TERUMO cannot guarantee the connection.

<Wireless LAN>

- 1) Use the TERUFUSION Drug Library Manager or TERUFUSION Software Package to set the communication environment in the pump.
- 2) Select a communication environment setting to be used for the pump.
- 3) Enable the wireless LAN of the pump.

Cautions


- When using the wireless LAN communication function, the effect on other equipment by radio interference should be considered.
- When connecting this product with a network, the settings of this product and the network system need to be compatible. Contact a vendor specialist for the correct settings and allow only the system administrator to perform settings.

Note

- TERUMO cannot guarantee the connection.

Transmission Procedure

When the network system sends an instruction (REQ: REQUEST) to the pump, the pump replies (REP: REPLY) to the instruction. The pump communicates with the network system on the predefined network (TCP Server) over the wireless LAN. The following shows an example of the communication sequence.

Notes

- It does not respond to data without a matching device ID (identification information), in other words, received content is not self-addressed.
- Wireless LAN communication is not available when it is attached to the rack system and is conducting IrDA communication.

About External Communication Functions

Packet

<Packet format>

The following shows the packet format.

Packet format

STX 0x02	Packet data (XML format) (Arbitrary size)	CRC16 (4byte)	ETX 0x03
-------------	--	------------------	-------------

(Packet data and CRC16 use ASCII code.)

- STX
0x02 (Fixed)
- Packet data
Data division written in XML format (For details, see this page "Packet <Packet data>".)
- CRC16
Error detection calculation for the packet data division. The polynomial to be used is CRC-16-CCITT.
- ETX
0x03 (Fixed)
- Measures for a packet error
If the pump receives data without STX and/or ETX, the data previously received become invalid.
If the pump receives data with the CRC error, the data previously received become invalid.

<Packet data>

For the XML declaration statement, only `<?xml version="1.0"?>` is allowed.

Do not use the DOCTYPE declaration.

A space, a tab, and any control code including a line break are not allowed.

XML entity reference notation is not supported.

• Request packet format

The following shows the request packet format used when a network system (or TERUFUSION Communication Rack System) makes a request to the pump.

Request packet format

<code><?xml version="1.0"?></code>			
<code><REQ></code>	<code><DEVID></DEVID></code>	<code><CMD></CMD></code>	<code></REQ></code>

- REQ tag:
This tag and its element are required.
This indicates a request packet from the network system to the pump.
The pump does not reply to a request without this tag.
- DEVID tag:
This tag is required.
This tag indicates the device ID of the pump.
The pump replies when the device ID of the element matches that of the pump.
When the network system obtains the device ID of the pump, the element is omitted. If element is omitted, the pump replies to only the DEVID tag without conditions, and sends its own device ID as the reply data.
- CMD tag:
This tag indicates the command and data to the pump.
The pump replies to one CMD tag.
For details on the element commands and data, see page 93 "Command Format <Request>".
The element is converted using Base64.

• **Reply packet format**

The following shows the reply packet format used when the pump replies to a request packet.

Reply packet format

<?xml version="1.0"?>			
<REP>	<DEVID></DEVID>	<CMD></CMD>	</REP>

- **REP tag:**
This tag indicates the reply from the pump to the request command sent by the network system.
- **DEVID tag:**
This tag indicates the device ID of the pump.
- **CMD tag:**
This tag indicates the command and data from the pump to the request command.
For details on the element commands and data, see this page "Command Format <Reply>".
The element is converted using Base64.

Command Format

<Request>

The following shows the format of a request command received by the pump.
The request command is used as the element of the CMD tag of a request packet format.

Request command format

Sequence	Command	SubCommand	Data
----------	---------	------------	------

- **Sequence:**
For a command from the network system: 0x00-0x7E
For a command from the TERUFUSION Communication Rack System: 0x80-0xFE
The following sequence numbers are reserved: 0xFF and 0x7F
- **Command:**
For details, see page 94 "Data Format".
- **SubCommand:**
For details, see page 94 "Data Format".
- **Data:**
Data is sent according to Command or SubCommand.
For details, see page 94 "Data Format".

<Reply>

The following shows the format of a reply command the pump sends.
The reply command is used as the element of the CMD tag of the reply packet format.

Reply command format

Sequence	Command	SubCommand	Execution result	Data
----------	---------	------------	------------------	------

- **Sequence:**
The pump sends the same value as the sequence number of the request command received.
- **Command:**
The pump sends the same value as that of Command of the request command received.
- **SubCommand:**
The pump sends the same value as that of SubCommand of the request command received.

About External Communication Functions

• Execution result:

When the pump has received a request command and successfully processed it, the pump sends "0x00" indicating that execution is successful.

If the request command failed due to an error, the pump sends a value between 0x11 and 0x15 as the execution result.

The following table lists the execution results.

Item	Size	Details
Execution result	1	0x00: Successfully processed
		0x11: Command or SubCommand cannot be identified or is undefined.
		0x12: The command is invalid, being unable to be executed
		0x13: Short message
		0x14: Long message
		0x15: The parameter is outside the range.

• Data:

If the execution is successful (0x00), the pump sends the processing results according to the requested command and the sub command.

If the execution failed due to the error (0x11), the pump sends 0x000000.

If the execution failed due to errors (0x12-0x15), the pump sends the requested command and the sub command.

Data Format

<Requesting basic operating status>

• Overview

The current basic operating status of the pump is obtained.

• Command : 0x01

• SubCommand : 0x0000: Basic operating status: Request

• Data

Network system → Pump : None

Pump → Network system : With data (See below)

The following table lists the details of the "Basic operating status" data:

Basic operating status (Do not handle reserved data)

Item	Size (Byte)	Details	
Software Ver.	2	Range: 0x0000-0xFFFF	
Basic operating status	1	Upper 4 bits Lower 4 bits	
		0x0	Stop (Reserved)
		0x1	Start (Reserved)
		0x2	Standby (Reserved)
		0x3	Priming (Reserved)
		0x4	Hands On Bolus (Reserved)
		0x5	Hands Free Bolus (Reserved)
		0x6-0xF	(Reserved) (Reserved)
(Reserved)	1	(Reserved)	
(Reserved)	1	(Reserved)	

Item	Size (Byte)	Details		
Alarm information	3	b0 (low order)	Shutdown Notice	
		b1	Link Interruption alarm	
		b2	Completion alarm	
		b3	Occlusion alarm	
		b4	Syringe Displacement alarm	
		b5	Plunger Displacement alarm	
		b6	Slider Displacement alarm	
		b7	Syringe Barrel Detection alarm	
		b8	Nearly Empty alarm	
		b9	(Reserved)	
		b10	(Reserved)	
		b11	(Reserved)	
		b12-b14	(Reserved)	
		b15	AC power detector buzzer sound function	
		b16	No Flow Rate alarm	
		b17	Flow Rate/VTBI Volume Judgment alarm	
		b18	No VTBI alarm	
		b19	Battery alarm	
		b20	Start Reminder	
		b21	Pressure alarm	
		b22	Maintenance timer	
		b23 (high order)	Battery Replacement Time/Failure alarm	
		One bit used for each item. ("1" indicates active.)		
		Power status	1	Upper 4 bits (Power type)
Lower 4 bits (battery level)				
0x0	Battery drive			Battery charge: 0 bars (All turned off)
0x1	AC drive			1 bar (less than 20%)
0x2	(Reserved)			2 bars (less than 40%)
0x3	(Reserved)			3 bars (less than 60%)
0x4	(Reserved)			4 bars (less than 80%)
0x5	(Reserved)			5 bars (greater than or equal to 80%)
0x6-0xF	(Reserved)	(Reserved)		
Operation monitor	1	b0-b1 (low order)	Occlusion monitor (0x1: 25%, 0x2: 50%, 0x3: 75%)	
		b2	(Reserved)	
		b3	Alarm sound status on/off	
		b4	Keypad lock status	
		b5	Plunger status monitor	
		b6	(Reserved)	
		b7 (high order)	(Reserved)	
		One bit used for b2 - b7. ("1" indicates active.)		

About External Communication Functions

Item	Size (Byte)	Details	
Syringe status	1	0x0	None
		0x1	(Reserved)
		0x2	5mL
		0x3	10mL
		0x4	20mL
		0x5	30mL
		0x6	50mL
		0x7-0xFF	(Reserved)
Flow rate	3	Value of 100 times the set flow rate	
Volume delivered	3	Value of 100 times the volume delivered	
VTBI	3	Value of 10 times the VTBI Infinite: 0xFFFFF Invalid value (when the function is disabled): 0xEEEEEE	
Dosage rate	3	Value of 100 times the dosage Invalid value (when the function is disabled): 0x000000	
Unit of dosage rate	1	0x00: mL/h 0x04: µg/kg/min 0x08: mg/kg/h	
Weight	2	Value of 10 times the weight Invalid value (when the function is disabled): 0x0000	
Dilution	3	Value of 100 times the dilution Invalid value (when the function is disabled): 0x000000	
Unit of dilution	1	0x13: mg/mL Invalid value (when the function is disabled): 0x00	
Drug library Ver.	2	Version of the drug library Drug library is registered: 0x0001-0xFFFF Drug library is not registered: 0x0000	
Drug library number	2	Number of the drug library Unit switching No.1 (µg/kg/min): 0x0FFD Unit switching No.2 (mg/kg/h): 0x0FFE Invalid value (when the function is disabled): 0x0000	
Occlusion pressure setting	1	bit0-bit3: Setting of occlusion pressure. 1 to 10. Occlusion pressure range: 0x01-0x0A bit7: Setting of occlusion detection pressure 1 = In 10 levels	
Hands On Bolus flow rate	2	100.00 to 1200.00 mL/h (100 mL/h step)	
Hands Free Bolus VTBI	3	Value of 100 times the bolus VTBI Invalid value (when the function is disabled): 0x000000	
Hands Free Bolus dose time (Bolus VTBI time)	2	0 to 3600 s (1 s step) Invalid value (when the function is disabled): 0x0000	
Hands Free Bolus flow rate	3	Value of 100 times the flow rate Invalid value (when the function is disabled): 0x000000	
Syringe brand number	1	Brand number of the syringe selected	
Total number of bytes	46		

* For more details, contact TERUMO trained service technicians.

Format of Transmitted and Received Data

When a network system communicates with the pump, the monitor obtains the device ID from the pump identifying itself. Then, the monitor starts to communicate with the pump.

The following shows the format of communication data.

1. Device ID request (network system → pump)

Data example

```
[STX]<?xml version="1.0"?><REQ><DEVID></DEVID></REQ>8411[ETX]
```

Note) [STX]: 0x02, [ETX]: 0x03

* In this case, CRC is "8411".

2. Device ID reply (pump → network system)

Data example

```
[STX]<?xml version="1.0"?><REP><DEVID>000100700000001201000001</DEVID></REP>D38C[ETX]
```

Note) [STX]: 0x02, [ETX]: 0x03

* This is an example of the reply where the device ID is "000100700000001201000001".
In this case, CRC is "D38C".

3. Basic operating status request (network system → pump)

Data example

```
[STX]<?xml version="1.0"?><REQ><DEVID>000100700000001201000001</DEVID><CMD>AQEAAA==</CMD></REQ>  
A9A4[ETX]
```

Note) [STX]: 0x02, [ETX]: 0x03

* This is an example of the basic operating status request where the device ID is "000100700000001201000001", the sequence number is "01", the command is "01", and the sub command is "0000".
In this case, CRC is "A9A4".

4. Basic operating status reply (pump → network system)

Data example

```
[STX]<?xml version="1.0"?><REP><DEVID>000100700000001201000001</DEVID><CMD>AQEAAAABAAAAAAAAAQBUA  
AAAnEAAAAAAD6AAAAAAAAAAAAAAAAAAAAAiQBkAAAAAAAAAAAA</CMD></REP>1E69[ETX]
```

Note) [STX]: 0x02, [ETX]: 0x03

* This is an example of the reply where the pump is stopped.
In this case, CRC is "1E69".

Procedure for Connecting to Nurse Call

The pump can output a nurse call alarm by using the RS-232C.

1) Connect the pins of the external communication/nurse call connector.

When alarm sounds (short)	6th pin, 9th pin
When alarm sounds (open)	1st pin, 4th pin

Cautions

- If connecting with a nurse call, always check in advance with the electrical contractor for compatibility, and connect with a system for which the connection function has been confirmed.
- If connecting with a buzzer or lamp instead of a nurse call, it should be used with a capacity of 12V DC, 1A or less.

Notes

- TERUMO cannot guarantee the connection.
- If an alarm sounds, this product sends signals to the nurse call for approx. 5 seconds.
- If an error is displayed, this product sends a continuous signal.

In order to ensure the safe use of the syringe pump, make sure it is cleaned and inspected after use and stored in an appropriate environment. If you find any damage or malfunction, immediately stop using the pump and contact TERUMO trained service technicians.

Cleaning

Cautions

- Clean the AC inlet, slider and clutch on a regular basis. Failure to conduct the regular cleaning may result in damage, failure or malfunction to/of this product.
- If any parts (clamp, clutch part etc.) of this product need to be removed for cleaning, contact TERUMO trained service technicians.
- When disinfecting, do not use a sterilizer. Use a gauze (damp with antiseptic solution) to wipe the product, and wipe off the antiseptic solution with a gauze (damp with cold/lukewarm water), and then thoroughly wipe off any moisture with a dry soft cloth.
- Before cleaning, always turn off the power and disconnect the AC power cable. Otherwise it may cause an electrical shock hazard. When AC power is supplied by the optional rack system (Model: TE-RS700, TE-RS800, TE-RS811), first remove the pump from the rack system, then conduct cleaning.
- If any adhered drug solution is present, the solution delivery or alarm detection may not operate correctly or may interfere with the normal operation of the clamp part and/or movable part of the pole clamp. If any drug solution is present, clean immediately with a cotton swab.
- Do not wipe with organic solvent such as alcohol or thinner, or with povidone-iodine. Using organic solvent or any antiseptic solution other than those permitted for use may result in damage to or failure of this product.
- Do not use an electric dryer to dry this product. This product may be damaged.
- This product does not have a waterproof structure. Therefore, do not wash with running water or immerse in water. It may cause damage or failure of the product.
- If the flange holder or clutch suffers from stiff movement or gets very dirty, it should be cleaned with a damp (not wet) soft cloth.

Examples of disinfectants (ingredient names) which can be used are listed below.

Examples of Cleaning Disinfectants (Ingredient Names)

Ingredient Name	Dilution (e.g.)
Chlorhexidine gluconate	5%
Benzalkonium chloride	10%

When using disinfectant, follow the instruction manual of each disinfectant (regarding degree of dilution, etc.)

How to Clean Components


After use, clean all components as described below.

Cleaning the AC inlet

Lightly wipe off dirt with a cotton swab or damp cloth (not wet).

Caution

- After cleaning, completely dry the AC inlet before use. Otherwise an electric shock or short circuit may occur resulting in failure of the product.


Cleaning the panel surface

Wipe with a damp (not wet) soft cloth.

Caution

- Do not rub the panel surface with a hard object.

Cleaning the outside of the syringe pump and the Dial


Wipe with a damp (not wet) soft cloth.


When disinfecting, use a designated disinfectant.

The Dial is fixed to the body of the pump with a magnet. Remove it following the procedure below before cleaning.

How to remove the Dial

- 1 Turn the groove of the Dial to the rear of the device.
- 2 Place a flathead screwdriver in the groove of the Dial and push it outwards.
- 3 Remove the Dial and wipe clean the body and the Dial.
Grease is applied to the shaft. Do not wipe it off.
- 4 When you have finished cleaning, fit the Dial back onto the body.
- 5 After fitting, turn on the power.
Check that setting for flow rate etc. can be performed.


Cleaning the pole clamp


Wipe with a damp (not wet) soft cloth.

Caution

- Do not soak the pole clamp in water to clean it.


• Pole clamp (screw type)

• TERUFUSION One Touch Pole Clamp (Optional accessory)

Cleaning the clamp

Wipe with a cotton swab or damp cloth (not wet).


Cleaning the flange slit and the Holding flange.

Wipe with a damp (not wet) soft cloth.

Cleaning the rubber protector and around the rubber protector

Wipe with a damp (not wet) soft cloth.

Storage

After inspecting the pump, store it in an environment which conforms to the following storage conditions

Temperature: -20°C to 45°C

Relative humidity: 10% to 95% (non-condensing)

Ensure that the battery is charged in preparation for it's next use.

Warnings

- Since this product does not have an airtight structure, it should not be used or stored in an active gas environment (including sterilizer gas), nebulizer-sprayed environment, high-humidity environment, etc. It should not be submerged into water. If the electronic components inside the product are affected, there may be subsequent damage and time degradation which will cause failure of this product.
- The product should not be used or stored in a flammable environment.

Cautions

- Do not store this product in a place with a high level of vibration, dust, mist or corrosive gas.
- Do not expose this product to direct sunlight or ultraviolet irradiation for a long period of time. The exterior may experience colour change, deformation or deterioration.
- Do not store this product in a place where an environment of atmospheric pressure, temperature, humidity, ventilation, salt content or sulfur content may cause an adverse effect.
- Do not store this product in a storage area for chemicals or in a place that generates gases.
- When transporting this product, avoid any shock, vibration, dust, high temperature and high humidity. Transport conditions: Ambient temperature -20 to 60°C, Relative humidity 10 to 95% RH (non-condensing)
- If the battery is not used for a long period after being charged, the amount of charge decreases due to self-discharge. Do not store the battery in a discharge status. Storing the battery left in a discharge status may cause degradation, resulting in making it unavailable in emergency situations.

Maintenance and Inspection Items by TERUMO Certified Service Technicians

Cautions

- If any fault is observed, immediately stop using the product and contact TERUMO trained service technicians.
- Do not disassemble, make alterations (including actions that interfere with the functionality or performance such as taping the LCD or a movable part) to, or repair this product. This may result in failure, damage or device performance degradation of this product.
- Do not use any replacement parts other than those specified. The original functionality or performance of this product may not be achieved.

Periodic Maintenance

Perform a periodic maintenance inspection to ensure safe operation and the longest possible life of the syringe pump. Contact your local distributor for details on the nature and frequency of the maintenance inspection.

From TERUMO

- When ordering repairs or maintenance inspections, if there is a possibility of infection, disinfect the product in advance.

Battery replacement

Periodically replaced parts

Over a period of time, the battery will gradually deteriorate. Please see chart below for replacement schedule.

Part name	Elapsed years	Replacement reason
Battery (internal battery)	2.5 to 3 years	Even when charged, after a short period of time "Battery alarm" is displayed and the buzzer sounds.

- * Depending on frequency of use and usage environment, the time for replacement for each part may change and replacement of parts other than those scheduled for replacement may also be required. Consult TERUMO trained service technicians regarding necessity and replacement.
- * The life span of the equipment is at least 6 years in the case of standard use.
- * This product incorporates a fuse inside the structure. For replacement, contact TERUMO trained service technicians.

Waste and Recycle

Electrical and electronic equipment (EEE) and batteries contain materials, components and substances which can be dangerous to the environment and harmful to human health if waste electrical and electronic equipment (WEEE) and batteries are not disposed of correctly.

Waste electrical and electronic equipment and batteries must not be disposed of with the remainder of unseparated waste, but should instead be collected separately. In this way, the environmental impact associated with disposal of WEEE and batteries is reduced and there will be more opportunity for reusing, recycling and recovering WEEE and recycling batteries.

The Li-ion battery should be removed from the equipment by TERUMO certified service technicians. Please contact your local distributor.

At end of life, please dispose of this equipment and batteries according to your local regulations. Contact your local distributor or municipality for details of the available collection schemes. The embedded memory back-up lithium coin battery, collected together with this equipment at end of life as WEEE, is to be removed and treated by the recycling center.

Troubleshooting

When an alarm/error occurs, the factors below may be the cause. Follow the procedure for dealing with each alarm/error as it arises.

Caution

- If the fault persists even after the following measures, a device failure is suspected. Contact TERUMO trained service technicians to request repair.


1) Alarm (e.g. High priority alarm).

Alarm displayed on LCD. Operation indicator flashes red.
Check the status such as stopped solution delivery.

Note


The priority of the alarm status is as follows:

- (1) High priority: Alarms flashing red on LCD
- (2) Medium priority: Alarms flashing yellow on LCD

2) Mute and check alarm content.

When the alarm is muted with the Back/Mute switch, the alarm content will be displayed. Press the Back/Mute switch again to return to the flow rate screen. (The type of alarm flashes in the lower part of the screen)


3) Close the three way stopcock etc. (if in use)

4) Inspect the cause. Check any applicable occurrence/cause from pages 105 to 116.


5) Eliminate the cause. Eliminate the cause in accordance with the instructions from pages 105 to 116.

6) Open the three way stopcock etc. (if in use)


7) Confirm and restart.

Check the flow rate, VTBI, etc
and resume the solution delivery.

◆ When installing a syringe: Plunger displaced


Place the syringe

Occurrence	Cause	Action
LCD: [Place the syringe.] Slider displacement icon flashes	The plunger has been displaced.	Place the syringe correctly. (See page 27 "Setting the Syringe".)

◆ Plunger Displacement alarm


⚠ Plunger Detection

⚠ Plunger Detection

Occurrence	Cause	Action
- LCD: Plunger Displacement alarm display flashes red - Operation indicator: Flashes red - Buzzer: Activates - Pump: Stops	The plunger has been displaced during solution delivery.	- Re-set the syringe correctly. - Prime the line, ensuring that the plunger is kept in the slider to prevent the solution from being delivered to the patient.
- LCD: Plunger Displacement alarm display flashes red - Operation indicator: Flashes red - Buzzer: Activates - Pump: Remains stopped, unable to start even after the syringe has been correctly set.	Device failure (plunger) A failure in the plunger detector is suspected. Failure is due to an external cause such as adhered drug solution or an internal cause such as sensor failure. After the alarm is issued, the solution delivery cannot start while the power is on.	- If any adhered drug solution is present, turn off the power and clean the plunger. - If the problem persists, contact TERUMO trained service technicians to request repair.

◆ When installing a syringe: Syringe displaced


Place the syringe

Occurrence	Cause	Action
LCD: [Place the syringe.] Syringe displacement icon  flashes	- The clamp is still open. - The flange is not in the slit. - The syringe is not installed correctly.	Place the syringe correctly. (See page 27 "Setting the Syringe".)
LCD: Syringe displacement icon  flashes	The syringe brand setting is incorrect.	Check that the syringe brand setting matches the brand of the syringe to be set.

◆ Syringe Displacement alarm


⚠Syringe Detection

⚠Syringe Detection

Occurrence	Cause	Action
- LCD: Syringe Displacement alarm display flashes red - Operation indicator: Flashes red - Buzzer: Activates - Pump: Stops	The syringe has been displaced during solution delivery.	- Replace the syringe correctly. - Prime the line, ensuring that the plunger is kept in the slider to prevent the solution from being delivered to the patient.
- LCD: Syringe Displacement alarm display flashes red - Operation indicator: Flashes red - Buzzer: Activates - Pump: Remains stopped, unable to start even after the syringe has been correctly set.	Device failure (flange or/and syringe clamp) A failure in the plunger detector or syringe detector is suspected. Failure is due to an external cause such as adhered drug solution or flange holder fault, or an internal cause such as sensor failure. After the alarm is issued, the solution delivery cannot start while the power is on.	- If any adhered drug solution is present, turn off the power and clean the flange or/and syringe clamp. (See page 100 "How to Clean Components".) - If the problem persists, contact TERUMO trained service technicians to request repair.

◆ When installing a syringe: Slider displacement

Place the syringe

Occurrence	Cause	Action
LCD: [Place the syringe.] Slider displacement icon flashes	The slider has been displaced.	Place the syringe correctly. (See page 27 "Setting the Syringe".)

◆ Slider Displacement alarm

⚠Slider Detection

Slider not in the correct position!
Clamp the line and check the slider position

⚠Slider Detection


Occurrence	Cause	Action
- LCD: Slider Displacement alarm display flashes red - Operation indicator: Flashes red - Buzzer: Activates - Pump: Stops	The slider has been displaced during solution delivery.	- Replace the syringe correctly. - Prime the line, ensuring that the plunger is kept in the slider to prevent the solution from being delivered to the patient.

◆ Occlusion alarm

Occurrence	Cause	Action
• LCD: Occlusion alarm display flashes red • Operation indicator: Flashes red • Buzzer: Activates • Pump: Stops	• An occlusion has occurred in the infusion line.	1. Eliminate the cause of occlusion. 2. Prime the line, ensuring that the plunger is kept in the slider to prevent the solution from being delivered to the patient.
• The occlusion icon does not turn red even when an Occlusion alarm is issued.	• The internal pressure of the infusion line has been temporarily mitigated by the bolus reduction function. • The cause of the occlusion may have been eliminated.	• Ensure that the cause of the occlusion is eliminated before resuming delivery. The bolus reduction function is not able to eliminate the cause of occlusion. Resuming delivery without eliminate the cause of occlusion may result in the solution delivery not being carried out correctly, including a repetitive Occlusion alarm or drug solution in the infusion line flowing back to the syringe.
• An Occlusion alarm is occurring often.	• The internal pressure of the infusion line has reached the occlusion pressure set value.	• If using a drug solution with a high viscosity or using a thin line, the internal pressure of the infusion line may increase. Review the settings by setting an occlusion pressure value higher than the one currently used, increasing the inner diameter of the infusion line, or reducing the flow rate in accordance with the occlusion icon or the occlusion history graph, etc.

◆ Pressure alarm (TE-SS800 only)

Occurrence	Cause	Action
• LCD: The Pressure alarm display flashes yellow • Operation indicator: Flashes green • Buzzer: Activates • Pump: Delivers solution	• The internal pressure of the infusion line has rapidly increased or decreased.	• Due to a possible fault in the infusion line or at the puncture site, check for any fault. Stop the solution delivery and eliminate the cause.


Note

- This alarm is effective only when the Pressure alarm has been set in the drug library.

◆ Nearly Empty alarm


⚠ Nearly Empty


Occurrence	Cause	Action
• LCD: Nearly Empty alarm display flashes yellow • Operation indicator: Flashes green (during solution delivery) or off (stopped). • Buzzer: Activates • Pump: No change to the status	• The quantity of the drug solution is low.	• If you want to continue the solution delivery, replace the syringe. • For replacing the syringe, see page 41 “Replacing Syringe to Continue Solution Delivery”.

◆ Syringe Empty


⚠ Syringe Empty

Occurrence	Cause	Action																		
• LCD: Syringe Empty display flashes red • Operation indicator: Flashes red • Buzzer: Activates • Pump: Stops	• The syringe is fully pressed.	• Replace the syringe if you wish to continue the infusion. The time from full-press to alarm sounding varies according to the product used. The table below shows an example. • When new TERUMO syringes are used. <table><tr><th></th><th>5 mL syringe</th><th>10 mL syringe</th><th>20 mL syringe</th><th>30 mL syringe</th><th>50 mL syringe</th></tr><tr><td>1 mL/h</td><td>Approx. 3 min.</td><td>Approx. 11 min.</td><td>Approx. 30 min.</td><td>Approx. 41 min.</td><td>Approx. 64 min.</td></tr><tr><td>5 mL/h</td><td>Approx. 1.0 min.</td><td>Approx. 2.5 min.</td><td>Approx. 5.0 min.</td><td>Approx. 8.5 min.</td><td>Approx. 12.0 min.</td></tr></table>		5 mL syringe	10 mL syringe	20 mL syringe	30 mL syringe	50 mL syringe	1 mL/h	Approx. 3 min.	Approx. 11 min.	Approx. 30 min.	Approx. 41 min.	Approx. 64 min.	5 mL/h	Approx. 1.0 min.	Approx. 2.5 min.	Approx. 5.0 min.	Approx. 8.5 min.	Approx. 12.0 min.
	5 mL syringe	10 mL syringe	20 mL syringe	30 mL syringe	50 mL syringe															
1 mL/h	Approx. 3 min.	Approx. 11 min.	Approx. 30 min.	Approx. 41 min.	Approx. 64 min.															
5 mL/h	Approx. 1.0 min.	Approx. 2.5 min.	Approx. 5.0 min.	Approx. 8.5 min.	Approx. 12.0 min.															


◆ No Flow Rate alarm

⚠No Flow Rate

Occurrence	Cause	Action
• LCD: No Flow Rate alarm display flashes yellow • Operation indicator: Off • Buzzer: Activates • Pump: Remains stopped	• The Start switch was pressed without having set a flow rate.	• Set a flow rate and press Start switch to start solution delivery.
	• The flow rate setting exceeds 150 mL/h (for 5 mL syringe) or 300 mL/h (for 10, 20 or 30 mL syringe).	• Check the value set and re-set it to equal or less then 150 mL/h or 300 mL/h, depending on the syringe size used.

◆ No VTBI alarm



⚠No VTBI

Occurrence	Cause	Action
• LCD: The No VTBI alarm display flashes yellow • Operation indicator: Off • Buzzer: Activates • Pump: Remains stopped	• The Start switch was pressed without having set a VTBI.	• Input a VTBI setting and press the Start switch to start solution delivery.
	• The Start switch was pressed with the settings as volume delivered $\geq$ VTBI.	• Clear the volume delivered or input a new value in the VTBI and press the Start switch to start solution delivery.

* Only if the VTBI setting function is enabled.


◆ Flow Rate / VTBI Volume Judgment alarm

Occurrence	Cause	Action
• LCD: The Flow Rate / VTBI Volume Judgment alarm display flashes yellow • Operation indicator: Off • Buzzer: Activates • Pump: Remains stopped	• The Start switch was pressed with the settings as Flow rate ≥ VTBI.	• Check the setting values of the flow rate and VTBI (if used). If they are not correct, input correct values and start solution delivery. If they are correct, press and hold the Start switch to start.

* Only if the VTBI setting function is enabled.



◆ Completion alarm

Occurrence	Cause	Action
• LCD: The Completion alarm display flashes yellow • Operation indicator: Flashes green • Buzzer: Activates • Pump: Switches to the keep vein open function and continues solution delivery.	• The volume delivered has reached the VTBI.	• A normal operation when the VTBI is set. • Press the Stop switch to stop the solution delivery.



* Only if the VTBI setting function is enabled.

◆ Syringe Barrel Detection alarm

Occurrence	Cause	Action
- LCD: Syringe Barrel Detection alarm display flashes yellow - Operation indicator: Off - Buzzer: Activates - Pump: Stopped	- The syringe in use does not match the syringe brand set. - The clamp is still open. - The flange is not in the slit. - The syringe has not been set at the center of the syringe detector.	- Check that the syringe brand setting matches the syringe in use. - Place the syringe correctly. (See page 27 "Setting the Syringe".)

◆ Start Reminder

Occurrence	Cause	Action
- LCD: The Start Reminder display flashes yellow - Operation indicator: Off - Buzzer: Activates - Pump: Remains stopped	The product has been left for 2 minutes or more after being programmed.	- When you do not start operation for a while, press and hold the Back/Mute switch for 2 seconds or more to set the equipment to standby mode. (See page 40) - Check the setting values for flow rate and VTBI, and press the Start switch to start solution delivery.

◆ Battery alarm

⚠ Low Battery

Low battery!
Connect the pump to the mains to
avoid shut down within 30
minutes

⚠ Low Battery

Occurrence	Cause	Action
- LCD: Battery alarm display flashes yellow. - Operation indicator: Flashes green (during solution delivery) or off (stopped). - Buzzer: Activates. - Pump: Continues solution delivery (during solution delivery) or remains stopped (stopped).	Low battery capacity (approximately 30 minutes), battery level low or battery energy is reducing.	- Connect to the AC power supply or rack system. * When using a rack system, see the instruction manual of the rack system.

◆ Shutdown Notice

⚠ Shut Down

Occurrence	Cause	Action
- LCD: Shutdown Notice display flashes red - Operation indicator: Flashes red - Buzzer: Activates - Pump: Stops	Battery has no power left.	- Connect to the AC power supply or rack system * When using a rack system, see the instruction manual of the rack system.

◆ Re-alarm

Occurrence	Cause	Action
- LCD: Respective alarm displays flash (the flashing colours are the same as those at the respective alarms). - Operation indicator: Same as for when the alarms are first activated - Buzzer: Activates - Pump: No change to the status	Operation was not started within 2 minutes after the alarm was muted.	Follow the corrective action of each alarm.

◆ Maintenance timer

Occurrence	Cause	Action
• The Maintenance timer indicator  lights up.	• It is time to perform maintenance.	• Contact TERUMO trained service technicians.

◆ AC power detector function

Occurrence	Cause	Action
• The pump cannot be turned on when using AC power.	The AC icon  is not displayed even after the AC power cable is connected.	• Check the connection of the AC power cable to the pump and the power supply.
	• The AC power cable is damaged.	• Replace the AC power cable.

◆ Battery Replacement Time/Failure alarm

Occurrence	Cause	Action
• The battery operation time is short even after fully charged. • The battery icon  is displayed. • The battery icon remaining battery life display number lights up to only 4 bars or less (after charging for 8 hours or more).	• The battery has deteriorated or is damaged.	• Replace the internal battery.

◆ Flow rate

Occurrence	Cause	Action
• Pressing the Start switch does not start solution delivery.	• The infusion line is not open. • The syringe is not installed correctly.	1. Replace the syringe correctly. 2. Prime the line, ensuring that the plunger is kept in the slider to prevent the solution from being delivered to the patient.
	• An occlusion has occurred in the infusion line.	• Eliminate the cause of occlusion (See page 108 "Occlusion alarm").
• Flow rate accuracy is low.	• Syringe has not been installed correctly.	• Please check the setting of the syringe. (See page 27 "Setting the Syringe".)
	• The syringe in use does not match the syringe brand set on the pump.	• Check that the syringe brand setting matches the brand of the syringe in use.
	• The flow rate is not set correctly.	• Correct the flow rate setting.
	• There is large height difference between the pump and patient.	• Adjust the patient to the pump height differential.

◆ External communication

Occurrence	Cause	Action
Cannot communicate (RS-232C)	The baud rate, stop bit, parity, or device ID number may not be appropriate.	Check the settings of the connected device. (See page 89 "About External Communication Functions".)
	Cables are not securely connected.	Reconnect the cable. (See page 89 "About External Communication Functions".)
Cannot communicate (IrDA)	The baud rate, stop bit, parity, or device ID number may not be appropriate.	Check the settings of the connected device.
	The infrared communication window.	Clean the infrared communication window with cotton swabs or gauze, etc.
	The pump is not securely attached to the TERUFUSION Communication Rack System.	Attach the pump properly.
	The TERUFUSION Communication Rack System is not turned on.	Turn on the TERUFUSION Communication Rack System.
Cannot communicate (wireless LAN)	SSID, encryption setting, IP address, subnet mask, default GW, connection destination IP address, communication port, or device ID may not be suitable.	Check the settings of the product or the connected device.
Signal strength icon (Colour = Grey, No signal bar) flashes	The pump is defective.	Contact TERUMO trained service technicians.

◆ Display

Occurrence	Cause	Action
The LCD is dim.	- The equipment is operated using the internal battery - The equipment is in the night mode. (When the night mode is enabled)	- The brightness becomes two levels lower during internal battery operation or night mode operation. - If it is not the case, contact TERUMO trained service technicians.
Operation indicator: Flashes red and green alternately.	The product is in standby mode.	Press the Start switch to start and standby mode will be cancelled.
The keypad lock icon  is displayed.	Keypad lock is on.	Keypad lock can be cancelled by holding down the Back/Mute switch (for 2 seconds or more).
Displays are irregular and the buzzer sounds.	The pump is defective.	Contact TERUMO trained service technicians to request repair.

◆ Failure

Occurrence	Cause	Action
"Er*" is displayed, the operation indicator lights up red and the buzzer sounds. (* is alphanumeric characters.)	The pump is defective.	Contact TERUMO trained service technicians to request repair.

Error Code Table

Display	Error	Display	Error
E1	Internal electronic circuit fault	E5*	Internal electronic circuit fault
E2*	Internal electronic circuit fault	E8*	Operation key fault
E3	Motor fault	E10*	Internal electronic circuit fault
E4*	Internal electronic circuit fault		

(* is alphanumeric characters.)

Caution

- Before pressing the Start switch to resume solution delivery after cancelling an alarm condition, ensure that you check the flow rate and VTBI settings.

Note

- If there is an alarm condition as the power is turned on, the alarm display flashes but the buzzer doesn't sound.

Main unit

Product name	TERUFUSION Syringe Pump Type SS
Model	TE-SS700/TE-SS800
Catalogue number	TE-SS700xxx, TE-SS702xxx, TE-SS800xxx (xxx = EN1, EN2, DE1, FR1, IT1, IT3, NL1, RU1, ES1, SV1, NO1, DA4, PL1, HU1, SK1) (See page 122 "Units Available" for the catalogue number.)
Compatible syringes	Sterilized specified syringes for general use and specified pre-filled syringes (See page 18 "Compatible Syringes")
Syringe size	5mL, 10mL, 20mL, 30mL, 50/60mL
Syringe brand	TERUMO or other specified brands
Dose mode (*: TE-SS800 only)	mL/h mode µg/kg/min mode (*) mg/kg/h mode (*) Library mode (*)
Flow rate setting range	- Setting range 0.01 to 150.00 mL/h Note that the upper limit of the flow rate can be changed in the following ranges: 0.01 to 150.00 mL/h (when using syringe of 5 mL) 0.01 to 300.00 mL/h (when using syringe of 10, 20, 30 mL) 0.01 to 1200.00 mL/h (when using syringe of 50/60 mL) - Step 0.01 mL/h step (0.01 to 10.00 mL/h)* 0.10 mL/h step (10.00 to 100.00 mL/h)* 1.00 mL/h step (100.00 to 1200.00 mL/h)* *: Step can be used when flow rate is set by the Dial. (It will be 0.00 mL/h when the power is turned on.)
VTBI setting range	0.10 to 100.00 mL (0.10 mL step) 100.00 to 9999.00 mL (1.00 mL step) When the dose rate is set in mass units, setting in mass units is possible. (VTBI setting range: 0.01 ng to 9999.99 kg) "----" (Unlimited) (It will be 0.00 mL when the power is turned on.)
VTBI time setting range	1 min to 99 hours 59 mins (1 min step) "--" (Unlimited)
Dose rate setting range (TE-SS800 only)	0.01 to 10.00 (0.01 step) 10.00 to 100.00 (0.10 step) 100.00 to 999.00 (1.00 step) (In units such as µg/kg/min and mg/kg/h etc.) (It will be 0.00 when the power is turned on.)
Weight setting range (TE-SS800 only)	0.1 to 300.0kg (0.1kg step)
Dilution setting range (TE-SS800 only)	0.01 to 10.00 (0.01 step) 10.00 to 100.00 (0.10 step) 100.00 to 999.00 (1.00 step) (In units such as mg/mL etc.)
Hands On Bolus flow rate setting range	100.00 mL/h (when using syringe of 5mL) 100.00 to 300.00 mL/h (when using syringe of 10mL, 20mL, 30mL) 100.00 to 1200.00 mL/h (when using syringe of 50/60mL) (100.00 mL/h step)

Specifications

Hands Free Bolus dose rate setting range		0.01 to 10.00 (0.01 step) 10.00 to 100.00 (0.10 step) 100.00 to 999.00 (1.00 step) (In units such as mL and mg/kg etc.)
Hands Free Bolus dose time setting range		1 s to 60 mins 00 s (1 s step)
Volume delivered display range		0.00 to 10.00 mL (0.01 mL step) 10.00 to 100.00 mL (0.10 mL step) 100.00 to 9999.00 mL (1.00 mL step) When the dose rate is set in mass units, the volume is expressed in mass units. (Volume delivered display range: 0.01 ng to 9999.99 kg) (It will be 0.00 mL when the power is turned on.)
Flow rate accuracy	Machine accuracy	Within $\pm 1\%$ (flow rate ≥ 1.00 mL/h)
	Accuracy including syringe	Within $\pm 3\%$ (The hourly precision after one hour has passed since the start of infusion at a constant flow rate of 1.00 mL/h or more according to IEC 60601-2-24 standard.)
Occlusion detection pressure		10 to 120kPa (set value) The occlusion detection pressure can be set in 10 levels.
Purge flow rate		Approx. 150 mL/h (when using syringe of 5mL) Approx. 300 mL/h (when using syringe of 10mL) Approx. 400 mL/h (when using syringe of 20mL) Approx. 500 mL/h (when using syringe of 30mL) Approx. 1200 mL/h (when using syringe of 50/60mL)
Alarms		Occlusion alarm, Pressure alarm* ¹ , Nearly Empty alarm, Slider Displacement alarm, Syringe Barrel Detection alarm, Syringe Displacement alarm, Plunger Displacement alarm, Battery alarm, Re-alarm, Start Reminder, No Flow Rate alarm, No VTBI alarm* ² , Flow Rate/VTBI Volume Judgment alarm* ² , Completion alarm* ² , Link Interruption alarm* ¹ * ¹ : TE-SS800 only * ² : Only if the VTBI setting function is enabled (Default: disabled)
Safety functions	Remaining volume detection setting function	Sets the Nearly Empty alarm position based on the time before the syringe is empty.
	Occlusion detection pressure select function	Selects the occlusion detection pressure setting values in 10 levels.
	Bolus reduction function	When an occlusion occurs, the slider is drawn back to reduce internal pressure.
	Flow rate range setting function A (Soft limits)	The flow rate setting range can be controlled within the "flow rate setting range" prescribed for each type of syringe. • When you try to set a value outside the range controlled by this function, a message asking for confirmation is displayed. If you confirm the setting, you can set the value which is outside the controlled range and the display colour of the value changes.
	Flow rate range setting function B (Hard limits)	The flow rate setting range can be controlled within the "flow rate setting range" prescribed for each type of a syringe. • Values outside the range controlled by this function cannot be set.
	Keypad lock function	Disables switch operations. However, the Start switch, the keypad lock function disable operation (hold down the Back/Mute switch), and the Power switch still function. When an alarm is activated, the keypad lock is released.

Available Functions	Syringe installation guidance function	Displays a screen that indicates the syringe installation status on the LCD.
	Syringe brand display function	Displays the syringe brand and size on the LCD.
	Syringe type switch function	Switching syringes is possible by performing the specified operation after the power is turned on.
	VTBI setting function	Sets the VTBI of dosage of drug solution.
	VTBI time setting function	Sets the VTBI time of dosage of drug solution.
	Purge/bolus buzzer sound function	When the Purge switch is being pressed, or during a bolus, the buzzer sounds.
	Volume delivered clear function	Clears the volume delivered display to zero.
	Standby function	By pressing the Back/Mute switch, the pump can be put in standby mode. In the standby function, the Start Reminder is delayed.
	Standby duration time select function	Changes the time to activate a Start Reminder in the standby mode.
	Buzzer volume select function	Changes the buzzer volume.
	Switch operation buzzer sound function	Buzzer sounds when a switch operation is performed.
	Stop transition buzzer sound function	Buzzer sounds when the Stop switch is pressed during solution delivery.
	LCD brightness select function	Changes the LCD brightness.
	Maintenance timer function	Time to next maintenance check can be set. That information is displayed on the LCD.
	Date and time setting function	Sets the date and time.
	External communication function (RS-232C) (TE-SS702 only)	Allows communication with a network system to output the status of the pump (flow rate, alarms,etc.) via RS-232C.
	Free message function (TE-SS702 and 800 only)	Displays a message on the pump via communication.
	External communication function (infrared communication (IrDA)) (TE-SS800 only)	Allows communication with the optional TERUFUSION Communication Rack System to output the status of the pump (flow rate, alarms,etc.) via infrared.
	External communication function (wireless LAN) (TE-SS800 only)	Allows communication with a network system to output the status of the pump (flow rate, alarms,etc.) via wireless LAN.
	Nurse call function (TE-SS702 only)	When an alarm is activated, relay contact operates through communication unit.
	History function	Displays the operation history.
	AC power detector function	Detects the presence or absence of AC power supply to the pump and displays that status.
	AC power detector buzzer sound function	Buzzer sounds when AC power supply is provided to or disconnected from the pump.

Specifications

Available Functions	Keep Vein Open Function	After activation of the Infusion Completion alarm, solution delivery is continued at the keep vein open flow rate.
	Keep Vein Open flow rate select function	Changes the keep vein open flow rate.
	During-start volume delivered clear function	Clears the volume delivered to zero during solution delivery.
	During-start flow rate/dosage rate change function	Changes the flow rate/dosage rate during solution delivery.
	Alarm melody select function	Changes the buzzing pattern (melody) of the alarm buzzer.
	Setting step 100 times function	When the Dial is turned while the Stop switch is pressed, setting steps can be performed at a rate of 100 times. (If the flow rate range setting function A is enabled and when you set the flow rate/dosage rate, this function is available only outside the flow rate range specified in the flow rate range setting function A.)
	Hands On Bolus function	After the bolus flow rate is set, a bolus is delivered while the Purge switch is pressed.
	Hands Free Bolus function	When bolus volume and bolus time are set, the flow rate is calculated and a bolus is delivered.
	Shortcut call function	Press multiple keys in combination to directly call the bolus flow rate setting display.
	Night mode function	Changes the LCD brightness for a specified time period.
	Power OFF buzzer sound function	Buzzer sounds as the power is turned off.
	Dose mode select function (TE-SS800 only)	The dose mode can be selected from the mL/h mode, µg/kg/min mode, mg/kg/h mode, and drug library mode.
	Patient information display function (TE-SS800 only)	Patient information including patient name, age, gender, weight, height, etc., is displayed.
	Interval/Intermittent dose function (TE-SS800 only)	The specified dosage rate is delivered at the specified regular intervals.
	Multi-step dose function (TE-SS800 only)	The speed of solution delivery increases or decreases at the specified levels (up to 10 levels).
	Delayed Start dose function (TE-SS800 only)	Solution delivery starts after the specified time has elapsed.
	Link function (Switch function A, B) (TE-SS800 only)	Switch function A: Operations of two pumps are linked. When a Completion alarm occurs on the first pump, the second pump starts solution delivery. Switch function B: Interlocks two pumps. Starts solution delivery of the first pump, and then starts solution delivery of the second pump at the predetermined time. After solution delivery of the second pump is started, each pump changes the solution delivery speed according to respective settings. The interlock of the pumps finishes upon completion of the set solution delivery or at alarm timing.
	Come and See ME function (TE-SS800 only)	Automatic call function after a specified time elapses.

Available Functions	Volume delivered unit select function (TE-SS800 only)	When the dose rate is set in mass units, the volume delivered unit can be changed to mass units.
	Micro notation select function (TE-SS800 only)	Switches "µg" to "mcg" to be displayed in the drug library mode, µg/kg/min mode or mg/kg/h mode.
Use conditions		Ambient temperature: 5 to 40°C Relative humidity: 20 to 90%RH (non-condensing)
Storage conditions		Ambient temperature: -20 to 45°C Relative humidity: 10 to 95%RH (non-condensing)
Transport conditions		Ambient temperature -20 to 60°C Relative humidity 10 to 95%RH (non-condensing)
Power supply		AC100 to 240V, 50 to 60Hz Internal battery (lithium ion battery) - Continuous use period: Approx. 12 hours (applies when solution is delivered continuously at 5 mL/h in ambient temperature of 25°C with new fully charged battery) - Charging time: ≥8 hours (when charged with AC power supply with the power turned off)
Power consumption		24VA
Classification		Class I equipment and internally powered equipment, type CF applied part, Continuous Operation, IP24
Wireless LAN		Hereby, TERUMO CORPORATION declares that this TERUFUSION Syringe Pump Type SS is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. - Transmit frequency or frequency band: 2412 – 2472 MHz - Type of modulation and frequency characteristics: - Differential Binary Phase Shift Keying: DBPSK - Differential Quadrature Phase Shift Keying: DQPSK - Direct Sequence Spread Spectrum: DSSS - Complementary code keying: CCK - Effective radiated power: 10.7 dBm
Dimensions		Without external communication (RS-232C): 381 mm (W)*3 x 120 mm (H)*3 x 112 mm (D)*4 With external communication (RS-232C): 381 mm (W)*3 x 120 mm (H)*3 x 131 mm (D)*4 *3: Excluding protrusions, *4: Excluding protrusions, pole clamp and moving range
Unit weight		Approx. 2.0 kg
Standard accessories		AC power cable, pole clamp (Screw type), instruction manual

* This product is in compliance with EMC (electromagnetic compatibility) standard IEC 60601-1-2:2001 Amd. 1:2004 (CISPR group classification and class classification are Group 1 and Class B). It is also in compliance with the EMC level required by IEC 60601-2-24:1998.

* Conformity standard
IEC 60601-1:1988, Amd. 1:1991, Amd. 2:1995
IEC 60601-1-1:2000
IEC 60601-1-2:2001, Amd. 1:2004
IEC 60601-1-6:2004
IEC 60601-1-8:2003, Amd. 1:2006
IEC 60601-2-24:1998
MDD (Medical Device Directive) 93/42/EEC (Class IIb)

Specifications

Units Available

Model	Catalogue number	Labeling	Plug type	Computer interface	Weight mode, Drug library mode
TE-SS700	TE-SS700EN1	English	European plug	-	-
	TE-SS700EN2	English	UK plug		
	TE-SS700DE1	German	European plug		
	TE-SS700FR1	French			
	TE-SS700IT1	Italian	European plug		
	TE-SS700IT3		Italian plug		
	TE-SS700NL1	Dutch	European plug		
	TE-SS700RU1	Russian			
	TE-SS700ES1	Spanish			
	TE-SS700SV1	Swedish			
	TE-SS700NO1	Norwegian			
	TE-SS700DA4	Danish	Danish plug		
	TE-SS700PL1	Polish	European plug		
	TE-SS700HU1	Hungarian			
	TE-SS700SK1	Slovak			
	TE-SS702EN1	English	European plug	SS702: yes (RS-232C)	-
	TE-SS702EN2	English	UK plug		
	TE-SS702DE1	German	European plug		
	TE-SS702FR1	French			
	TE-SS702IT1	Italian	European plug		
	TE-SS702IT3		Italian plug		
	TE-SS702NL1	Dutch	European plug		
	TE-SS702RU1	Russian			
	TE-SS702ES1	Spanish			
	TE-SS702SV1	Swedish			
	TE-SS702NO1	Norwegian			
	TE-SS702DA4	Danish	Danish plug		
	TE-SS702PL1	Polish	European plug		
	TE-SS702HU1	Hungarian			
	TE-SS702SK1	Slovak			
TE-SS800	TE-SS800EN1	English	European plug	SS800: yes (IrDA, Wireless LAN)	SS800: yes
	TE-SS800EN2	English	UK plug		
	TE-SS800DE1	German	European plug		
	TE-SS800FR1	French			
	TE-SS800IT1	Italian	European plug		
	TE-SS800IT3		Italian plug		
	TE-SS800NL1	Dutch	European plug		
	TE-SS800RU1	Russian			
	TE-SS800ES1	Spanish			
	TE-SS800SV1	Swedish			
	TE-SS800NO1	Norwegian			
	TE-SS800DA4	Danish	Danish plug		
	TE-SS800PL1	Polish	European plug		
	TE-SS800HU1	Hungarian			
	TE-SS800SK1	Slovak			

Note

- Separate setup procedures are required for the product with voltage-plug combinations other than listed above.

Test data for flow rate characteristics and occlusion characteristics is shown below.
Tests were conducted based on IEC 60601-2-24: 1998. For details, refer to IEC 60601-2-24: 1998.


Flow Rate Characteristics

Start-up curves and trumpet curves, which are used as typical explanations of flow rate characteristics, are shown. The data below is representative examples measured for flow rates of 1 mL/h and 5 mL/h with a new TERUMO syringe (50mL).


Note

- Flow rate characteristics are correlated with the syringe. Depending on the combination of syringe and pump, the same data as the representative example may not be obtained.

- Start-up curve
The discharge quantity measured every 30 seconds for a 2-hour measurement period is converted into flow rate values is shown. These graphs show characteristics from immediately after the start of solution delivery to when flow rate has stabilized.

- Trumpet curve
Data for the latter one hour of the 2-hour measurement period is shown, with the horizontal axis showing the observation window time (minutes) and the vertical axis showing the maximum and minimum flow rate error for each observation window. The narrower the region enclosed between the top and bottom solid lines (the so-called trumpet curve) is, the less pulsatory motion there is.

Occlusion Characteristics

Typical data of occlusion pressure detection ability are occlusion pressure, occlusion alarm occurrence time, and post-occlusion removal bolus volume.

The data below is a representative example measured for flow rates of 1 mL/h and 5 mL/h, with a three-stage occlusion pressure setting value (1, 6, 10), and a new TERUMO syringe (50 mL).

* When the 1m extension tube (internal volume: 1.0 mL) is used.

Flow rate (mL/h)	Occlusion pressure setting value	Occlusion pressure		Time to activate an Occlusion alarm time (min)	Bolus volume (mL)
		(kPa)	(kgf/cm ²)		
1	1	13	0.14	18.1	0.05
	6	70	0.72	71.2	0.34
	10	119	1.21	134.6	0.83
5	1	13	0.14	2.3	0.07
	6	72	0.73	14.7	0.33
	10	120	1.22	22.2	0.81

EMC Technical Information

Medical electrical equipment requires particular care in regards to EMC, and it is necessary to install and use the equipment in accordance with the following EMC information.

Cautions

- Ensure that the power cable used is the one provided for this product. If not, this product's emissions may increase and its immunity may decrease.
- When using this product in combination with other equipment such as a network system, do not place the equipment next to each other or stack them together. It may cause malfunctioning due to electromagnetic interference.

Table 201 – Guidance and manufacturer's declaration – electromagnetic emissions –
for all EQUIPMENT and SYSTEMS (see 6.8.3.201 a) 3))

Guidance and manufacturer's declaration – electromagnetic emissions		
The TERUFUSION Syringe Pump Type SS is intended for use in the electromagnetic environment specified below. The customer or the user of the TERUFUSION Syringe Pump Type SS should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The TERUFUSION Syringe Pump Type SS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class [B]	The TERUFUSION Syringe Pump Type SS is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class [A]	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 202 – Guidance and manufacturer's declaration – electromagnetic immunity –
for all EQUIPMENT and SYSTEMS (see 6.8.3.201 a) 6))

Guidance and manufacturer's declaration – electromagnetic immunity			
The TERUFUSION Syringe Pump Type SS is intended for use in the electromagnetic environment specified below. The customer or the user of the TERUFUSION Syringe Pump Type SS should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±8kV contact ±15kV air (*)	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TERUFUSION Syringe Pump Type SS requires continued operation during power mains interruptions, it is recommended that the TERUFUSION Syringe Pump Type SS be powered from an uninterruptible power supply or a battery.
Power frequency (50-60 Hz) magnetic field IEC 61000-4-8	3 A/m	400 A/m (*)	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note

- U_T is the a.c. mains voltage prior to application of the test levels.

* The TERUFUSION Syringe Pump Type SS complies with the more stringent levels of IEC 60601-2-24 : 1998.

Table 203 – Guidance and manufacturer’s declaration – electromagnetic immunity –
for LIFE-SUPPORTING EQUIPMENT and SYSTEMS (see 6.8.3.201 b))

Guidance and manufacturer’s declaration – electromagnetic immunity			
The TERUFUSION Syringe Pump Type SS is intended for use in the electromagnetic environment specified below. The customer or the user of the TERUFUSION Syringe Pump Type SS should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands ^a	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the TERUFUSION Syringe Pump Type SS, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P}$
	10 Vrms 150 kHz to 80 MHz in ISM bands ^a	10 Vrms	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	$d=1.2 \sqrt{P}$ 80 MHz to 800 MHz $d=2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^c , should be less than the compliance level in each frequency range ^d . Interference may occur in the vicinity of equipment marked with the following symbol: 
Notes 1. At 80 MHz and 800MHz, the higher frequency range applies. 2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a. The ISM (industrial, scientific and medical) bands between 150 KHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27. 283 MHz; and 40.66 MHz to 40.70 MHz. b. The compliance levels in the ISM frequency bands between 150 KHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges. c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TERUFUSION Syringe Pump Type SS is used exceeds the applicable RF compliance level above, the TERUFUSION Syringe Pump Type SS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the TERUFUSION Syringe Pump Type SS. d. Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 205 –Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM–for LIFESUPPORTING EQUIPMENT and SYSTEMS (see 6.8.3.201 b))

Recommended separation distances between portable and mobile RF communications equipment and the TERUFUSION Syringe Pump Type SS				
The TERUFUSION Syringe Pump Type SS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TERUFUSION Syringe Pump Type SS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TERUFUSION Syringe Pump Type SS as recommended below, according to the maximum output power of the communications equipment.				
Separation distance according to frequency of transmitter				
m				
Rated maximum output power of transmitter W	150 kHz to 80 MHz outside ISM bands $d=1.2 \sqrt{P}$	150 kHz to 80 MHz in ISM bands $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.3
100	12	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
Notes 1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. 2. The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. 3. An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. 4. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

This section describes the symbols (icons) used with the TERUFUSION Syringe Pump Type SS.

Description of Symbols

Symbol	Description of Symbols	Symbol	Description of Symbols	Symbol	Description of Symbols
	TERUFUSION Syringe Pump Type SS is in conformance with the Medical Device Directive, 93/42/EEC.		Start		Relative humidity limitation
	TERUFUSION Syringe Pump Type SS is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.		Stop		Stacking limit by number
	Authorized Representative in the European Community		IEC 60529 Degrees of protection provided by enclosures (IP code)		Type CF applied part
	Serial number		Indicates the connecting position of External communication/nurse call (RS-232C)		Indicates recyclable corrugated paper packing material (Japan)
	Catalogue number		Indicates the operation to pull the clamp forward		Indicates PE (polyethylene) plastic bag (Japan)
	Date of manufacture		Contents		Recyclable (lithium-ion battery) Brazil
	Manufacturer		Consult instructions for use		Recyclable (lithium-ion battery) China
	Alternating current		This way up		Recyclable (lithium-ion battery) Taiwan
	"off" (for a part of equipment)		Fragile, handle with care		Recyclable (lithium-ion battery) Japan
	"on" (for a part of equipment)		Keep away from rain		Separate collection of accumulators and batteries
	Return		Keep away from sunlight		Separate collection of electrical and electronic equipment, European Community
	Mute (Audio paused)		Temperature limitation		Non-ionizing electromagnetic radiation
	Purge				

Note

- See pages 5 to 11 for symbols on the operation panel and LCD.

FOR INFORMATION ABOUT TERUMO PRODUCTS

- If this product should fail to perform as intended, immediately stop using the product and contact the nearest branch or sales office of TERUMO.

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TERUMO DENMARK FILIAL OF TERUMO SWEDEN AB

Tel.: +45 7020 93 80
Fax: +45 7020 94 80

TERUMO UK Ltd.

Tamesis, The Causeway
Egham, Surrey, TW20 9AW
United Kingdom
Tel.: +44 1784 476 200
Fax: +44 1784 476 201

TERUMO EUROPE N.V. MOSCOW REPRESENTATIVE OFFICE

BC Northern Tower
13Fr., 10 Testovskaya Street
Moscow 123317, Russia
Tel.: + 7 495 988 4740
Fax: + 7 495 988 4739

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® : Registered trademark ©TERUMO CORPORATION 2013-01-01 13A18