FotoFinder

Videodermatoscope

"FotoFinder bodystudio ATBM"

User Manual

Please read this manual carefully before use and
always keep it easily accessible!

2016

Page 1 of 122

FotoFinder Systems GmbH

Industriestraße 12
D-84364 Bad Birnbach
Germany

Tel.: +49 (0) 8563 - 97720-45

Fax: +49 (0) 8563 - 97720-10

Manufacturer authorized representative

MAXmedical, LLC

6 Sirenevaia Street, Office 54
Moskovskaia oblast, 141109 Shchjol'kovo
Russian Federation

Phone: +7 499 502 79 44,

Fax: +7 499 502 50 36.

www.maxmedical.ru

E-mail:

General questions: info@maxmedical.ru

Tech. support: support@maxmedical.ru

Content

Part I. Technical manual.

1. Videodermatoscope "FotoFinder bodystudio ATBM"	7
1.1. Components of the videodermatoscope	8
2. General and technical information	9
2.1. Universally used descriptions	9
2.2. Conventional signs	12
2.3. Intended purpose of the Videodermatoscope "FotoFinder bodystudio ATBM"	13
2.4. Conditions of use	14
2.5. Users	15
2.6. Proper handling rules	15
2.6.1. Putting into service	15
2.6.2. Installation instructions	16
2.6.3. Environmental conditions	16
2.6.4. Information on repair	17
2.7. Cleaning and disinfection	17
2.8. Servicing and safety	18
2.8.1. Videodermatoscope "FotoFinder bodystudio ATBM"	18
2.8.2. Leaser Liner USB	19
2.8.3. Medicam 800HD	19
2.9. Handling of the photo and video cameras	19
2.9.1. Photo cameras Canon EOS 700D, EOS 750D	19
2.9.2. Medicam 800HD	20
3. Installation at a working place	21
3.1. Preparation for operation	21
3.1.1. Photo cameras Canon EOS 700D and EOS 750D settings	21
3.1.2. Connecting photo cameras Canon 700D and Canon 750D to PolFlash face device and FotoFinder Server	21
3.1.3. Preparations for photography	22

4. Technical data	22
4.1. Videodermatoscope "FotoFinder bodystudio ATBM"	22
4.2. FotoFinder Silent Medical Server	23
4.3. Laser Liner USB	24
4.3.1. Working with the Laser Liner USB	24
4.3.2. Marking	25
4.4. Medicam 800HD	26
4.4.1. Working with Medicam 800HD	26
4.4.2. Spacer	27
4.4.3. D-Scope II Plus Lens.....	28
4.4.4. Medicam 800HD controls	30
4.4.5. FD-Lens	31
4.4.6. D-Scope III Lens.....	32
5. Working in the computer network	32
6. Troubleshooting	33
7. Testing	36
7.1. Technical and safety check-up	36
7.2. Visual inspection	36
7.3. Operational lifetime of the components of the videodermatoscope	36

Part II. Software manual.

1. Videodermatoscope "FotoFinder bodystudio ATBM"	38
2. Starting the system and opening the FotoFinder universe program	39
3. Section "Patients"	40
3.1. Adding new patient	40
3.2. Amending patient records	42
3.2.1. Opening a patient record	42
3.2.2. Making patient data anonymous	42
3.2.3. Editing patient information	43
3.2.4. Deleting a patient record	43

4. Section "Total Body Mapping" (TBM)	44
4.1. The interface	45
4.2. Photography	50
4.2.1. Taking baseline pictures	50
4.2.2. Additional body segments	52
4.2.3. Taking follow-up pictures	53
4.3. Image treatment and analysis	56
4.3.1. Editing image	57
4.3.2. Comparing images	59
4.3.3. Navigational tools	60
4.3.4. Image preview tools	63
4.3.5. Bodyscan ATBM	66
4.3.5.1. In the Compare mode	66
4.3.5.2. When adding markers	69
4.3.5.3. On the TBM desktop	69
4.3.6. Inserting and copying a marker	70
4.3.7. Connection to Dermoscopy	74
4.3.7.1. Adding a dermoscopy image	74
4.3.7.2. Displaying dermoscopy images	76
4.3.7.3. Comparing dermoscopy images	79
5. Section "Dermoscopy"	80
5.1. User interface	80
5.1.1. Main menu	81
5.1.2. Patient data area	82
5.1.3. Image preview window	82
5.1.4. SmartGallery	83
5.2. Overview pictures	86
5.2.1. Overview pictures with photo camera Canon EOS	86
5.2.2. Overview pictures with Medicam 800HD	87
5.2.3. Overview pictures with mannequin	89
5.2.4. Inserting marker	90
5.3. Dermoscopy images with Medicam 800HD	92
5.4. Images comparison	96
5.4.1. Comparing dermoscopy images	96

5.4.2. Comparing overview images	98
5.4.3. Other functions in the image comparison mode	100
5.5. Printing and reports	102
5.6. Optional: Moleanalyzer program	103
5.7. Optional: PASIvision program	105
5.8. Optional: FotoFinder Trichoscale Pro program	106
6. Section "Screening"	107
6.1. User interface	107
6.2. Main menu	107
6.3. Screening procedure	108
6.4. Saving images	108
7. Section "Query"	109
7.1. Full text search	111
7.2. Extended search	112
7.3. Sorter function	115
8. Settings	116
8.1. General	116
8.2. Imaging devices	117
8.3. Other devices	118
8.4. Modules	118
8.5. Changing user information	120
9. Security of data	121

Part I. Technical manual.

1. Videodermatoscope "FotoFinder bodystudio ATBM"

Automated Total Body Mapping with the FotoFinder bodystudio ATBM allows rapid, standardized documentation of the entire surface of the skin of patients at risk of developing skin cancer. The software guides the user through 20 specified positions. The system is further optimized to document Psoriasis. This technical manual gives you details on proper use, safety, maintenance and trouble shooting. Regarding the software refer to Part II "Software manual" on Page 38.


Photo 1. Videodermatoscope "FotoFinder bodystudio ATBM"

1.1. Components of the videodermatoscope


Photo 2. Components of the videodermatoscope.

2. General and Technical Information

2.1. Universally used descriptions

Original	Casually used descriptions
General	
FotoFinder Systems GmbH	<i>Full:</i> FotoFinder Systems GmbH <i>Short:</i> FotoFinder or: Manufacturer
Principal components of the videodermatoscope	
Videodermatoscope "FotoFinder bodystudio ATBM"	<i>Full:</i> Videodermatoscope "FotoFinder bodystudio ATBM". <i>Short:</i> Videodermatoscope.
FotoFinder Silent Medical Server	<i>Full:</i> FotoFinder Silent Medical Server. <i>Short:</i> FotoFinder Server.
Medicam 800HD	<i>Full:</i> Medicam 800HD. <i>Short:</i> Medicam or: video camera
LCD Monitor	<i>Full:</i> LCD Monitor. <i>Short:</i> Monitor.
Isolating Transformer	<i>Full:</i> Isolating Transformer.
ATBM Tower Station	<i>Full:</i> ATBM Tower Station. <i>Short:</i> Tower Station.
Linear Drive "dryLin"	<i>Full:</i> Linear Drive "dryLin". <i>Short:</i> Linear Drive.
Photo camera Canon EOS without lens, manufacturing modifications EOS 700D, EOS 750D	<i>Full:</i> Photo camera Canon EOS without lens, manufacturing modifications EOS 700D, EOS 750D <i>Short:</i> Photo camera Canon EOS or: Canon EOS or casual: photo camera.
SIGMA 30 мм F1.4 DC HSM Lens	<i>Full:</i> SIGMA 30 мм F1.4 DC HSM Lens. <i>Short:</i> SIGMA Lens.

Laser Liner USB	<i>Full:</i> Laser Liner USB. <i>Short:</i> Laser Liner.
PolFlash face	<i>Full:</i> PolFlash face.
AC/DC Power Adapter with power cable	<i>Full:</i> AC/DC Power Adapter with power cable. <i>Short:</i> AC/DC Power Adapter or: Power Adapter.
Flash-Coupler	<i>Full:</i> Flash-Coupler.
Stepper motor	<i>Full:</i> Stepper motor.
Control Unit TSA	<i>Full:</i> Control Unit TSA.
Power Distribution Unit	<i>Full:</i> Power Distribution Unit.
Medical Network-Isolation	<i>Full:</i> Medical Network-Isolation.
Power cable 150 cm long	<i>Full:</i> Power cable 150 cm long. <i>Short:</i> Power cable.
Computer keyboard *	<i>Full:</i> Computer keyboard. <i>Short:</i> Keyboard.
Computer mouse *	<i>Full:</i> Computer mouse. <i>Short:</i> mouse.

* sourced locally

Accessories

Canon EF-S 18-55 mm Lens	<i>Full:</i> Canon EF-S 18-55 mm Lens. <i>Short:</i> Canon EF-S Lens.
D-Scope II Plus Lens	<i>Full:</i> D-Scope II Plus Lens. <i>Casual short:</i> the lens.
Polarizing Filter for D-Scope II Plus Lens	<i>Full:</i> Polarizing Filter for D-Scope II Plus Lens. <i>Casual short:</i> the polarizing filter.
Open top piece for D-Scope II Plus Lens	<i>Full:</i> Open top piece for D-Scope II Plus Lens. <i>Short:</i> Open top piece.


Closed top piece for D-Scope II Plus Lens	<i>Full:</i> Closed top piece for D-Scope II Plus Lens. <i>Short:</i> Closed top piece.
Open glass tube for D-Scope II Plus Lens	<i>Full:</i> Open glass tube for D-Scope II Plus Lens. <i>Short:</i> Open glass tube.
Closed glass tube for D-Scope II Plus Lens	<i>Full:</i> Closed glass tube for D-Scope II Plus Lens. <i>Short:</i> Closed glass tube.
Conical glass tube for D-Scope II Plus Lens	<i>Full:</i> Conical glass tube for D-Scope II Plus Lens. <i>Short:</i> Conical glass tube.
D-Scope III Lens	<i>Full:</i> D-Scope III Lens. <i>Casual short:</i> the lens.
FD-Lens	<i>Full:</i> FD-Lens <i>Casual short:</i> the lens.
Black cloth for FD-Lens	<i>Full:</i> Black cloth for FD-Lens <i>Short:</i> Black cloth.
FotoFinder Pol Body	<i>Full:</i> FotoFinder Pol Body. <i>Short:</i> Pol Body.
Positioning mat	<i>Full:</i> Positioning mat.
ATBM Photo Background (blue)	<i>Full:</i> ATBM Photo Background (blue). <i>Short:</i> Photo Background.

Programs and functional tools

FotoFinder universe program	<i>Full:</i> FotoFinder universe program. <i>Short:</i> FotoFinder universe or casual: the program.
Moleanalyzer program	<i>Full:</i> Moleanalyzer program. <i>Short:</i> Moleanalyzer or casual: the program.
FotoFinder Trichoscale Pro program	<i>Full:</i> FotoFinder Trichoscale Pro program. <i>Short:</i> Trichoscale Pro or casual: the program.

PASlvision program	<i>Full:</i> PASlvision program. <i>Short:</i> PASlvision or <i>casual:</i> the program.
Bodyscan ATBM module	<i>Full:</i> Bodyscan ATBM module. <i>Short:</i> Bodyscan ATBM or: Bodyscan.
Cam-control function	<i>Full:</i> Cam-control function. <i>Short:</i> Cam-control.
Smart Gallery module	<i>Full:</i> Smart Gallery module. <i>Short:</i> Smart Gallery.
Ghost Function	<i>Full:</i> Ghost function.
SetCustomerInfo instrument	<i>Full:</i> SetCustomerInfo instrument. <i>Short:</i> SetCustomerInfo.
Other definitions	
Dermoscopy image(s)	<i>Short:</i> micro image(s)

2.2. Conventional signs

Equipotential earthing

Connection to protective conductor

Observe the user instructions

Alternating current

General warning sign

Only suitable for the interior

Beware of strip intake

Symbol "Do not dispose of as home waste"

CE symbol

Symbol "Manufacturer"

Symbol "Laser danger"

2.3. Intended purpose of the Videodermatoscope "FotoFinder bodystudio ATBM"

Videodermatoscope "FotoFinder bodystudio ATBM" is intended to be used for automated photography of the entire skin surface of patients followed by analysis of the saved images for the purpose of locating skin lesions (mole mapping) and conducting their computerized dermoscopy analysis. The device also enables the documenting of skin surface by psoriasis estimating the Psoriasis Area and Severity Index (PASI) and enables imaging of hair covered spots on the head followed by the image analysis and by drawing-up a trichogram.

The videodermatoscope enables making the following procedures:

- Automated photo documentation of the entire surface of patients skin;
- Documenting of data related to patients' pictures;
- Skin lesions mapping and classification allowing to locate them on the body pictures and trace their changes;
- Non-invasive digital dermoscopy of intact patients' skin;
- Documenting and analysis of the affected area and severity of psoriasis;
- Making of pictures of the haired areas on the head and creating of trichogram.

The videodermatoscope may then be efficiently used for the following:

- Repeated documenting and long-term monitoring of the entire body skin surface targeting the detection of new and changed skin lesions;
- Fast screening for skin cancer;
- Computer dermoscopy analysis of dermoscopy images of skin lesions by Moleanalyzer program;
- Fluorescent diagnosis for the differentiation of precancerous and cancerous stages of non-melanocytic lesions from a healthy skin;
- Documenting of psoriasis and calculating of Psoriasis Area and Severity Index (PASI);
- Creating the digital trichogram by FotoFinder Trichoscale Pro application.

There are no known counter-indications to the use of the videodermatoscope.

2.4. Conditions of use

The videodermatoscope should be used only externally (non-invasive) for taking of automated photo images of the entire body skin surface and for the computer aided analysis of the images using mole mapping and dermoscopy sections of the FotoFinder universe program.

The videodermatoscope is intended to be used exclusively with the FotoFinder universe program.

Only the following components supplied by the manufacturer should be used with the videodermatoscope: Photo camera Canon EOS 700D or EOS 750D, SIGMA 30 mm F1.4 DC HSWM Lens and Canon EF-S 18-55 mm Lens. The said photo cameras and the lenses are not equipped with a polarizing filter. It is also recommended to use the ATBM Photo Background supplied by FotoFinder as a background for photography.

For any questions concerning optimal photography conditions please consult with the authorized representative of FotoFinder Systems in Russian Federation, MAXmedical, LLC at: 6 Sirenevaia Street, office 54, 141109 Moskovskaia oblast, Schyolkovo. Phone +7 499 502 79 44, fax +7 499 502 50 36. E-mail: support@maxmedical.ru .

General conditions of using the videodermatoscope:

- The distance between patient and the background should be about 15 to 20 cm. The distance to the background should be the same at follow ups.

- The distance from the ATBM Tower Station post to the red line on the Positioning mat should be about 113 cm. More precise aligning of the Positioning mat can be done against the red beam of the working Lase Liner.
- Only use the device in bright and well illuminated rooms. Avoid direct solar radiation.
- It is only possible to make pictures of patients with a height of 130-200 cm. Complete photography of larger and smaller patients is not supported.
- The videodermatoscope is designed and equipped for an average 8 hour long working day and the patient number of max. 20 patients per working day.
- Analysis with Bodyscan module can only be valid for adult patients (more than 18 years old), as it is otherwise not possible to guarantee correct documentation because of body height changes.
- The classification of the Bodyscan results is based on statistical analyses and does not replace a professional diagnosis. The diagnosis is an all-time responsibility of the doctor.

2.5. Users

Working with the videodermatoscope is only permitted for doctors. Several technical procedures may be referred to assistants trained by a doctor.

2.6. Proper handling rules

2.6.1. Putting into service

Installation and pre-starting procedure of the videodermatoscope "FotoFinder bodystudio ATBM" should be only made by the authorized representative of the manufacturer, FotoFinder Systems GmbH, MAXmedical, LLC.

Before the installation of the videodermatoscope, all questions of its correct placing and connection should be agreed upon with the customer, i.e.:

- Minimal size of the place for making the photography – 1,3 x 3,2 m
- Availability of the operable grounded outlet that should be not more than 1,0 m from where the videodermatoscope is to be installed.
- Availability of proper lighting necessary for providing standard working conditions in the room.

The following should be checked every time before switching the videodermatoscope on:

- Check the reliability of all cable connections and fix them if necessary;
- All components of the videodermatoscope should be connected to the equipotential earthing bus. The grounding system in the working room should function reliably and ensure protection against electric shock.

It is required to stop working on the videodermatoscope before checking up and fixing the malfunctions if:

- Any visible damages to cables have been discovered;
- Photo camera Canon EOS or Medicam 800HD have been dropped on a hard floor or experienced a strong shock or lack integrity;
- Electric or electronic components have visible damages or produce unusual noise or there is smell of burning or smoke coming out of them. In this case, it is necessary to immediately unplug the entire system from the power source!

2.6.2. Installation instructions

The videodermatoscope must be installed at its working place so that it could stand sturdy. The cart rollers have brakes to prevent accidental movement of the device.

It is important to watch that no foreign objects could get into the Linear Drive and its belt transmission. All cables (for example, from Medicam) should be put aside from the Linear Drive and properly fastened to prevent them getting accidentally in touch with the moving mechanisms.

Do not place any objects in the close vicinity of the Linear Drive to avoid its failure or a breakdown of the photo camera Canon EOS and mechanisms.

During the installation of the videodermatoscope it should be ensured the FotoFinder Server has an adequate air inflow (not less than 50 cm on any side) – this will help to prevent failures due to overheating.

Videodermatoscope should not be installed closer than 50 cm to heating batteries or 1,0 m to infrared heaters or in places where it could be subject to continuous direct sunlight (more than half an hour in summer or two hours in winter), excessive dust, strong vibrations or shocks.

2.6.3. Environmental conditions

The videodermatoscope is intended to be used in the medical and scientific medical institutions including hospitals, clinics, outpatient clinics and educational medical institutions. The use is allowed only inside buildings under the following conditions:

- Ambient temperature +10°C to +35°C;
- Relative humidity up to 80% under the temperature of +25°C;
- Atmospheric pressure from 84 to 106,7.

2.6.4. Information on repair

Making any changes to the configuration of the videodermatoscope by its user is not allowed. Upon discovering of any such changes the warranty will be immediately void!

The manufacturer does not assume any responsibility for any damages or malfunctions if the videodermatoscope was used incorrectly or if its configuration was changed in any way.

Never open the case of the Laser Liner because you may be subject to uncontrolled effect of the laser beam.

Before cleaning and checking of the Medicam it must be completely disconnected from the power source.

For all questions concerning malfunctioning or failure of the equipment please contact the authorized representative of FotoFinder Systems, MAXmedical, LLC at: Phone +7 499 713 35 55, fax +7 499 501 50 36. E-mail: support@maxmedical.ru.

2.7. Cleaning and disinfection

Cleaning of the videodermatoscope:

- Before cleaning of the system make sure it is completely disconnected from the power supply source.
- Clean the body, the panels and the control elements with a soft cloth lightly moistened with a non-aggressive cleaning remedy. Do not use any abrasive materials or coarse cleaning appliances and solvents like alcohol or petrol.
- Clean the LCD Monitor with a soft cloth. If you use glass cleaning stuff make sure it does not contain antistatic substances and other ingredients with the same properties because they may scratch the protective monitor coating.

Cleaning of the Medicam 800HD:

- Clean the body of the Medicam with a soft cloth lightly moistened with a non-aggressive cleaning remedy. Do not use any abrasive materials or coarse cleaning appliances and solvents like alcohol or petrol.
- The glass tubes for the D-Scope II Plus Lens should be disinfected before every use. If you use an alcohol based liquid as the immersion fluid the glass tube is getting disinfected automatically. However, that does not happen while using immersion oil. Therefore it is recommended to use alcohol based liquids. If you use the polarizing kit and work without immersion fluids, disinfect regularly the plastic top piece.

- To clean the glass tube please pull it out from the lens. Use quality glass cleaning stuff and a lint-free cloth to clean the inside and the outside of the tube. Wipe the tube thoroughly with a dry cloth afterwards.
- To clean the D-Scope II Plus Lens please moisten the cloth with a small amount of the glass cleaning stuff.
- If there is dust on the lens's glass get rid of it with a soft brush.
- It is recommended once in four to six months to get the rubber ring that holds the glass tube lubricated with a small quantity of mineral oil but please avoid spreading it onto the lens. While putting the glass tube into the lens please make sure it is inserted in a straight position and all the way down. Otherwise the autofocus of the video camera while taking the dermoscopy images may not work correctly.

2.8. Servicing and safety

2.8.1. Videodermatoscope "FotoFinder bodystudio ATBM"

Automatic Linear Drive is made of components that do not require maintenance. They need just a simple periodical cleaning with a brush and/or vacuum cleaner (it is recommended to do once a year).

It is important to watch that the photo camera console and the belt transmission could move easily and without jerks. To ensure that, watch that any cables (for example, the video camera cable) are fastened so that they could not get in touch with the moving mechanisms. Make sure that no foreign objects could get into the moving parts of the Linear Drive.

Do not place any objects in the close vicinity of the Linear Drive to avoid its failure and breakdown of the photo camera and mechanisms.

Never touch the mechanisms of the Linear Drive with hand. Always try to observe the distance of not less than 30 cm from the moving parts (photo camera Canon EOS and belt transmission).

If the movement of the photo camera Canon EOS was stopped during its automatic travel (for example with hand or by pressing the emergency stop button in the program interface), all pictures taken for this particular patient pose will be deleted and the session for this pose has to be repeated by pressing the "Back" button in the program and starting the session again.

While breaks during work the photo camera Canon EOS may slowly move down under its own weight. In this case it is recommended to bring the photo camera to the lowest position on the Linear Drive where it will be protected against the outside mechanical risks with a special protective shield.

2.8.2. Laser Liner USB

It is necessary to check the distance which the Positioning mat is placed at against the videodermatoscope. It is specified on the drawing located on the right side of the Laser Liner and should be 113 ± 2 cm from the foot of the ATBM Tower Station post to the red line on the mat marked "Laser line". The ultimate mat aligning is made according to the laser beam of the working Laser Liner. The measurement tool is not supplied with the videodermatoscope, therefore the distance should be measured with a roulette.

While working with the turned on Laser Liner it is not permitted to allow the direct impact of the laser beam to the eye retina – i.e., your gaze should be directed at the laser beam always at an angle!.

2.8.3. Medicam 800HD

During the idle time, the Medicam should always be hooked onto the special holder located on the right side of the drawer desk of the ATBM Tower Station. Medicam should never experience any violent mechanical shocks and falls. If due to a mechanical shock or any different impact the Medicam stopped working, please check its cable and the consistency of its connection to the FotoFinder Server, and if it does not work then contact for maintenance the authorized representative of the manufacturer, FotoFinder Systems GmbH, MAXmedical, LLC company at: 6 Sirenevaia Street, office 54, 141109 Moskovskaia oblast, Schyolkowo, Russian Federation. Phone +7 499 713 35 55. E-mail: info@maxmedical.ru.

Medicam 800HD is subject to prophylactic maintenance once a year by the said authorized representative of the manufacturer.

2.9. Handling of the photo camera and video camera

2.9.1. Photo camera Canon EOS 700D or EOS 750D

- Do never disassemble the photo camera.
- Keep track of the working condition and reliable connection of the Canon EOS cables.
- Do not allow any liquids to leak into the photo camera.
- Do not touch the lens of the photo camera.
- If you removed any lenses, do not place them with lens facing down to avoid it to be scratched,
- Before cleaning the photo camera, disconnect the system completely from the power source.

- If the photo camera produce any usual noises, smoke or odor, please disconnect the entire system from the power supply immediately and contact the authorized representative of FotoFinder Systems, MAXmedical, LLC, at: Phone +7 499 713 35 55. E-mail: support@maxmedical.ru .

2.9.2. Medicam 800HD

Medicam 800HD is an electronic software-controlled video camera, developed to take microscopic and macroscopic pictures of the intact patient's skin, especially of pigment moles.

Few general instructions before using the Medicam:

- Do not deconstruct the video camera.
- Do not twist the video camera cable or step on it.
- Do not spill any liquids on the video camera.
- Do not touch the lens of the video camera.
- Be careful when taking off the D-Scope II Plus Lens and putting it on a glass table or similar surfaces because the glass lens may easily be scratched.
- Before cleaning the video camera, please completely disconnect the system from the power source.
- If the video camera makes any unusual noise, produces smoke or unusual odor, please disconnect the entire system immediately from the power source and contact the authorized representative of FotoFinder Systems GmbH, MAXmedical, LLC company at: 6 Sirenevaia Street, office 54, 141109 Moskovskaia oblast, Schyolkowo, Russian Federation. Phone +7 499 713 35 55. E-mail: info@maxmedical.ru .
- When inserting a glass tube into the D-Scope II Lens, make sure to pull it in correctly. Wrong handling may prevent the integrated automatic focus of the video camera from working correctly.
- We strongly recommend using alcohol spray as the immersion liquid. When using immersion oil the inside of the tube might get dirty. Make sure to remove any oil residues so that the video camera may focus correctly.

3. Installation at a working place

The videodermatoscope is delivered to the end user largely assembled. It is only needed to mount few components on the ATBM Tower Station – FotoFinder Server, LCD Monitor and the PolFlash face device with the photo camera Canon EOS, and connect all cables. This job is be done by the authorized representative of FotoFinder Systems, MAXmedical, LLC.

3.1. Preparation for operation

3.1.1. Photo cameras Canon EOS 700D and EOS 750D settings

As standard, the videodermatoscope is delivered with a Canon EOS 700D or EOS 750D photo camera and a Sigma 30 mm F1.4 DC HSM Lens. The photo camera has optimal factory settings the use with the videodermatoscope. Nevertheless, you should make sure your photo camera use the following settings:

- Camera mode set by the menu wheel: M (manual);
- Software settings:
 - o Auto power off: OFF;
 - o Auto lighting optimizer: OFF.
- Sigma lens: AF-mode (autofocus).

3.1.2. Connecting photo cameras Canon 700D and Canon 750D to the PolFlash face device and the FotoFinder Server

1. Insert the PolFlash face device's built-in adapter into the battery compartment of the photo camera.
2. Connect the photo camera to the PolFlash face using the screw connection.
3. Install the PolFlash face with the photo camera on the Linear Drive moving console with the screw provided in it.
4. Insert the Flash-Coupler into the flash slot on top of the photo camera.
5. Install the lens onto the photo camera.
6. Connect the power supply cable from the AC/DC Adapter to the PolFlash face device.
7. Connect the black-and-yellow USB cable to the photo camera.
8. Switch on the photo camera by pressing its ON/OFF button next to the menu wheel.

9. Check if the menu wheel of the photo camera is set to mode: M.
10. Installation of the Canon EOS driver starts automatically as soon as the photo camera is connected to the working FotoFinder Server via the USB cable.

3.1.3. Preparations for photography

1. Remove the remaining packing stuff from the Linear Drive mechanisms.
2. Switch on the photo camera.
3. Switch on the Laser Liner.
4. Align the Laser Liner beam with the red line on the Positioning mat.
5. Start the FotoFinder universe program.
6. Help the patient to take right position on the mat according to the software's instructions.

4. Technical data

4.1. Videodermatoscope "FotoFinder bodystudio ATBM"

Manufacturer:	FotoFinder Systems GmbH
Address:	Industriestrasse 12, D-84364 Bad Birnbach, Germany.
Name:	Videodermatoscope "FotoFinder bodystudio ATBM"
Power voltage/frequency:	AC 100-230 V / 47-63 Hz
Power consumption:	350 Watt
Protection class:	1
IP protection class:	IPX0
Transport and storage:	The device is shipped upright on a euro-pallet by a professional carrier company.

Disposal and environmental protection:

The device cannot be disposed of as domestic waste. Please disposal of the product in a professional and environmentally friendly way.

Compatible Canon photo cameras: Canon EOS 700D, EOS 750D

Note: The devices are manufactured in accordance with ISO 9001:2008.

4.2. FotoFinder Silent Medical Server

Type:	PC (personal computer)
Power supply:	Single phase alternative current from Isolating Transformer.
Power voltage:	100-240 V
Power frequency:	50-60 Hz
Nominal power consumption:	250 VA
Processor:	Intel®Core™ i5, 2,4 Ghz
RAM:	Not less than 8192 Mb (8Gb)
Video card:	Not less than 2 Gb.
Disc drive:	CD/DVD
Hard drive:	Not less than 1 Tb
System solid drive (SSD):	Not less than 128 Gb
Internet access:	For updates
Operational System:	Windows® 8.1 Professional (64-bit), Russian
USB ports:	Not less than 6 free ports
Dimensions (length, width, height):	435 x 395 x 160 mm
Weight:	10,3 kg
Grounding:	The FotoFinder server has grounding wire in the power cable. There is a contact on the body where the equipotential earthing bus is connected to.

to.

4.3. Laser Liner USB

Name:	Laser Liner USB
Laser Liner Wavelength:	650 nm (visible spectrum)
Maximum average radiant power:	Not more than 5,0 mW
Beam diameter at laser aperture:	0,6 mm
Emission duration:	≥ 0,25 seconds
Laser Class:	1
Fan angle:	20°
Working distance:	1,5 m/ 59 inch
Frequency drift:	0.25 nm/ °Celsius
Input voltage:	5.0 VDC

Caution:

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

4.3.1. Working with the Laser Liner USB

Turn on the Laser Liner by pressing the small ON/OFF button at the rear panel of the device. The ON/OFF button is marked with the following symbol:

Turn off the device by pressing the same button again.

Notes:

- The Laser Liner turns off automatically 7 minutes after power on.
- The Laser Liner will be in non-stop mode if you press the on/off button for more than 5 seconds!

However, a too long operation is not recommended.

4.3.2. Marking


The following labels are pasted onto the laser product:

Left side:


A label with warning sign and text.

Right side:


This label shows the distance from the post of the ATBM Tower Station to the red line on the Positioning mat.

A label with the name of the device, the manufacturer and its address.

Front:

A warning label

Back:

The On/Of button

4.4. Medicam 800HD

Type:	Full HD digital video camera
Model:	Medicam 800HD
Protection class:	1
Application part type:	B
IP protection class:	IPX0
Image pickup device:	1/3 Type CMOS
Number of effective pixels:	≈2,000,000 pixels
S/N ratio:	More than 50dB
Zoom:	20x optical
Electronic shutter:	½ - 1/10.000 sec. 21 steps
Power source:	FotoFinder Server
Power voltage:	5V / 12 V DC
Power consumption:	No more than 10W
Type of lighting:	White LED, 18 ea
Luminous intensity:	1000 lx
Spectrum and color temperature:	White daylight, 4200K
Light intensity control:	No
Working distance to object:	20 mm to 800 mm
Field of view size:	1,0 m x 1,0 m

4.4.1. Working with Medicam 800HD

Medicam is completely controlled by the FotoFinder universe software. This allows you to take identical images without adjusting the video camera settings. This makes work with FotoFinder universe extremely easy, fast and efficient. The Medicam is best equipped for a stable non-stop use during each working day.

The release button is integrated in the video camera handle. Press this button to take an image.

Caution: Please be careful with the video camera cable! Do not twist it or stand on it.

Floodlight illumination

Your video camera is equipped as standard with an integrated floodlight illumination.


Photo 4. Video camera Medicam 800HD with built-in floodlight illumination

The integrated floodlight illumination allows you to take excellent overview images with a distance of up to 80 cm. To guarantee perfect light conditions please use a recommended background. Please use the distance holder for standardized images without camera vibrations. On the backside of the video camera you can switch the illumination on and off.

4.4.2. Spacer

Medicam 800HD is equipped with a distance holder to take clinical images.

Thanks to the detent button it is possible to choose between two fixed distances. Keep the small detent button pressed down and pull the distance holder into the right position until it locks.

The two fixed distances allow a standardization of close-up images and facilitate to measure images.

Note: Use the distance holder together with the integrated floodlight illumination for images at a distance of 2 to 6 cm.


Photo 5. Medicam 800HD with the D-Scope II Plus Lens.

4.4.3. D-Scope II Plus Lens

The D-Scope II Plus Lens can be attached to the video camera via a bayonet clutch. Simply put the lens onto the contact area of the video camera and turn the clutch clockwise until it locks. The lens can easily be removed by turning it anti-clockwise.

A LED powered illumination is integrated into the lens. These LEDs are durable for a life time and do not need to be changed. The illumination is switched on by pushing the button on the back panel of the video camera. When the lens is installed on the video camera the floodlight is automatically disconnected.

The glass tube is held in the lens's depression by a rubber ring and can be removed for cleaning by gently twisting it clockwise and pulling out.




Photo 6. The D-Scope II Plus Lens.

Accessories for performing the immersion dermoscopy:

Closed glass tube for D-Scope II Plus Lens for standard immersion dermoscopy.

Open glass tube for D-Scope II Plus Lens for immersion dermoscopy on prominent lesions.


Conical glass tube for D-Scope II Plus Lens for immersion dermoscopy in places with limited access.

Accessories for performing the polarized dermoscopy:

1. Closed top piece for D-Scope II Plus Lens;
2. Open top piece for D-Scope II Plus Lens for polarized dermoscopy on prominent lesions.
3. Polarizing filter for D-Scope II Plus Lens.

4.4.4. Medicam 800HD controls

The green light diode shows that the video camera is turned on and ready for use. Please use the ring marked button for images taken with the D-Scope II Plus Lens, D-Scope III Lens or the FD -Lens.

The floodlight button is used for overview and clinical images. It turns off automatically after two minutes.

4.4.5. FD-Lens

FD-Lens for fluorescence diagnostics delineates actinic keratosis and basal cell carcinomas from the surrounding healthy tissue. To make use of it, a photo sensitizer is applied to the patient's skin. After an incubation time of 3 to 5 hours porphyrin is induced in tumor cells. Under the deep violet light of the FD-Lens the porphyrin molecules get excited and the tumor cells become visible and can be delineated from the healthy skin.

Thanks to the lens's chrome spacer and the specially designed cloth the skin is completely covered when the image is taken. You can toggle between white and violet light using buttons on the back panel of the video camera.

Photo 7. FD-Lens.

Caution! The lens contains light emitting diodes. Do not stare into the light.

Procedure

1. Please wait 3 to 5 hours until the photo sensitizer is induced in tumor cells.
2. Attach the FD-Lens to the video camera by inserting the pins on the right and left side of the lens to the holes in the video camera's contact ring and turn the bayonet

clutch of the lens clockwise until it is fixed firmly. By turning the lens anti-clockwise it can be pulled off the video camera.

3. Select "Medicam" as an imaging device and press the "Overview" button. **Do not choose the micro image capturing mode!**
4. Turn on the ring light button on the back panel of the video camera.
5. Place the front of the FD-Lens onto the skin.
6. Take a clinical color image and then mark the localization on the mannequin before saving it.
7. Please use the black cloth together with the lens in order to ensure optimal imaging conditions.
8. Now switch to fluorescence light by pressing the left button on the camera's back panel and take a fluorescent image. Then you will need to specify the localization again.

4.4.6. D-Scope III Lens


Photo 8. D-Scope III Lens.

D-Scope III Lens enables making dermoscopy images of small areas of skin with high magnification – up to 250X. That can be necessary for detailed examination of micro structure of building elements of skin – for example, in capillaroscopy. D-Scope III Lens is built to be used with video camera Medicam 800HD.

5. Working in the computer network

The videodermatoscope can work in a computer network with several FotoFinder universe clients. To ensure a galvanic isolation the Medical Network-Isolation device can be plugged in between the network cable and the "Network" socket of the FotoFinder Server.

It is also possible to adjust the FotoFinder universe software for use with the former version, FotoFinder R3 software.

Please contact the authorized representative of FotoFinder Systems, MAXmedical, LLC at: 6 Sirenevaia Street, office 54, 141109 Schyolkovo, Moskovskaya oblast. Phone: +7 499 502 79 44, fax +7 499 502 50 36. E-mail: support@maxmedical.ru .

6. Troubleshooting

This section contains information on troubleshooting. Please try to work through the following points step by step. If none of the steps solves the problem and the system does not start functioning correctly again, then please contact the authorized representative of FotoFinder Systems, MAXmedical, LLC at: 6 Sirenevaia Street, office 54, 141109 Schyolkovo, Moskovskaya oblast. Phone: +7 499 502 79 44, fax +7 499 502 50 36. E-mail: support@maxmedical.ru .

Note. In case of software problems the support may also be done remotely by the technical team of the manufacturer, FotoFinder Systems GmbH. If you wish to use it, please download and install the TeamViewer software from the following site: www.fotofinder.de/support . Then simply send your ID and password shown in the TeamViewer window at: support@fotofinder.de . The program is also necessary for updates, even if you do not experience any software problems.

A. Problems with the software

Case 1: Software error: The software is no longer functioning correctly

- Close the program, then restart the FotoFinder Server and run the program again.
- If the error persists or the software is still not functioning correctly, please contact the FotoFinder support team.

Case 2: The emergency stop in the FotoFinder universe program has been pressed

- Close the FotoFinder universe program and then run it again.

B. Problems with the hardware

Case 1: Photo camera Canon EOS is not reacting to commands

- Check the ON switch on the photo camera and the control light on the PolFlash face.
- Restart the photo camera by switching it off and on again
- Check the connection cable (CamControl cable), USB cable, power cable and built-in power adapter in the photo camera battery compartment.
- Check the photo camera settings: Auto lighting optimizer must be disabled and Auto power off switched off. The SIGMA lens should be set to the AF (Auto Focus) mode.

Case 2: PolFlash face – the work indicator does not light

- Check if the PolFlash face device is switched on.
- Then check the flash base socket – the plug should be completely inserted.
- Check the power connector and the supply of the PolFlash face for damages and correct connection.

Case 3: The stepper motor is not working or the software shows an error notification that the tower control could not be initialized

- Check the connection cable on the rear panel of the FotoFinder Server.
- Restart the FotoFinder Server and the FotoFinder universe software.
- Remove the FotoFinder Server and, if necessary, the printer from the bottom box and open the cover of the bottom box. Now check the control light and cable connections on the Control Unit TSA.
- Check the ON switch on the Control Unit TSA, switch the device off, wait approx. 20 seconds and switch it back on again.
- Check all fuses in the Control Unit TSA.
- Open the Device Manager in the Windows Control Panel and check whether the device "Silicon Laps CP210x USB to UART Bridge (COM3)" exists at the COM&LPT port.

Case 4: FotoFinder Server does not start

- Check the ON switch of the power module on the FotoFinder server rear panel.

- Check if the power cable is connected to the power module correctly.
- Check the ON switch on the front of the FotoFinder Server (the control light should be illuminated orange).
- Press the ON switch on the front (the control light should now illuminate blue)
- To sort out any problems with electricity, connect the videodermatoscope to a different outlet.

Case 5: LCD Monitor has a black screen

- Check the ON switch on the front and rear panels of the LCD Monitor.
- Check the power cable and all connections between the LCD Monitor and the FotoFinder Server.
- Check the signal (DVI, VGA or Display Port), press the "S" button on the front of the LCD Monitor several times until it is set correctly.
- To sort out any problems with electricity, connect the videodermatoscope to a different outlet.

Case 6: Medicam 800HD does not work

- Check the green control light on the rear of the Medicam.
- Check the connections on the rear panel of the FotoFinder Server.
- Shut down the FotoFinder Server, disconnect and reconnect the video camera cables and then restart the FotoFinder Server.

Case 7: Printer does not work

- Check the ON switch on the front and rear panels of the printer.
- Check the USB cable between the printer and the FotoFinder Server.

Case 8: Isolating transformer: no power supply to the system

- Check the green control light on the Isolating Transformer.
- Check or replace both fuses. Replace them with same type fuses only!
- Check all cables and their connections.
- Try a different power cable.

7. Testing

7.1. Technical and safety check-up

Please check the system for visible damages every time before you use the system. It is especially important to check if the power supply lines or the video camera are damaged.

Please check the system regularly or at least once every 24 months.

7.2. Visual inspection

Please have your device checked and inspected at regular intervals. The inspection should be done at least once every 12 months. As a minimum, the following things should be checked during a visual inspection:

- Protection against direct contact with moving parts.
- Protection measures in case of indirect contact must be in accord with the installation standards described.
- Preventive measures against over-current.
- Valid measurements of protection measures against over-current. Cross-section of conductors should correspond to power consumed.
- Circuit diagrams, descriptions of electrical circuits, operating instructions.
- Safety markings and labels to be intact.
- Status of grounding system.
- The condition of power and communication cables.
- Power light on the Control Unit TSA.
- Condition of cart rollers.

7.3. Operational lifetime of the components of the videodermatoscope

The components of your ATBM system have a certain service life. To ensure you are aware of the expected service life, you can find more information below. However, please note that this information is based on reference values and imprecise estimates. The actual service life of individual components can differ from these specifications.

The service life is given below in imaging cycles, i.e., complete photo documentation of patients. It is assumed that an imaging cycle consists of 20 images taken by the photo camera.

Service life of the individual components (estimated values):

– PolFlash face:	24.000 cycles.
– The belt of the Linear Drive and the Stepper motor:	100.000 cycles.
– Photo camera series Canon EOS:	5.000 cycles.
– Hard drive of the FotoFinder Server:	8.000 cycles.

To ensure perfect functioning of your ATBM system, your device should be subject to inspections and safety&technical checks periodically.

The device should work correctly and safely. Primarily, all electrical components should be inspected visually.

None of the components of the Linear Drive (the belt, the photo camera moving console and the rail) requires maintenance or greasing. These components should be visually inspected for wear (the belt and the rail) and stretching (the belt) periodically.

If these components are no longer working well they must be replaced professionally. If you have any problems with the system of automated photography, please contact the authorized representative of FotoFinder Systems, MAXmedical, LLC at: 6 Sirenevaia Street, office 54, 141109 Schyolkovo, Moskovskaya oblast. Phone: +7 499 502 79 44, fax +7 499 502 50 36. E-mail: support@maxmedical.ru.


Some components need to be checked more often than others and maintained or repaired as necessary.

Please observe the above specified regulations governing the safety&technical test and safe operation of all electric and moving parts and components.

Part II. Software manual.

1. Videodermatoscope "FotoFinder bodystudio ATBM"

FotoFinder bodystudio ATBM (Automated Total Body Mapping) provides an easy and reliable means of conducting total body skin examinations, especially for high risk melanoma patients. The FotoFinder universe program offers a standardized imaging procedure that guides you through 20 different positions. Various features guarantee for a capturing process that is as straight-forward as possible: while the video camera is software-controlled, the Laser Liner helps to ensure the distance and the Positioning mat gives clear indications of how to position the patient.

Distinctive features:


- Simplified imaging process.
- Bodyscan module to automatically locate changed and suspect lesions.
- Direct connection to dermoscopy.

2. Starting the system and opening the FotoFinder universe program

FotoFinder universe program (Version 2.0.33.2 of 27.01.2016) is started from your Windows desktop. In order to run the program, simply double-click the **FotoFinder.Universe** icon.

When the program is started for the first time, the *End User Licence Agreement* window will pop up. Please read this agreement carefully. You will be able to proceed once you have accepted the agreement at the bottom of the page. After checking the box on the left, this window will not appear again in the future.

You will now be taken to the FotoFinder universe dashboard, the platform from which to start for the conduction of body skin examinations. It consists of the following eight areas and sections, which are available depending on your program licensing:

- | | |
|------------------|------------------------------|
| A Settings | E Section Screening |
| B Patient Search | F Section Dermoscopy |
| C Patients | G Section Aesthetics* |
| D Query | H Section Total Body Mapping |

While using FotoFinder universe program you can log out of a section and log into another one without having to close or reset the program via *Disconnect*.

*Note. The Aesthetics section (letter G on the above screenshot) is not enabled in this particular configuration of the videodermatoscope "FotoFinder bodystudio ATBM".

The below screenshot shows that the program switched to different section:


3. Section "Patients"

In order to run a Total Body Mapping, the program requires you to first select an existing patient record or to add a new patient record to the database.

3.1. Adding new patient

To add a new patient record, hit the *Add patient* button in the upper left hand corner of the Patients section.

The following form will appear:

Please fill in the patient information. Blanks are provided to fill in mandatory and optional patient information.

Fields marked with an asterisk are mandatory!

The format of the patient record number is not set up as default and can be defined by the user individually. It can consist of letters and numbers, or a combination of both. Automatic assignment of the patient record number can be selected in the *Settings (A)* under *General* by ticking the second box ("generate new patient record numbers automatically").

While body height is not marked as a mandatory field, this piece of information is indeed required to run the Total Body Mapping (TBM). If you do not fill in the field at this point, you will be asked to provide the body height information before the first exposure in the Total Body Mapping module.

When you have filled out the patient information form, click on **Save**. A new patient record has been added to the database and you are referred back to the dashboard.


Save

3.2. Amending patient records

3.2.1. Opening a patient record


In order to select an existing patient record for update, you have two possibilities. You can either search the record in the search mask on the Dashboard itself, or log into the *Patients* section.

In either case, insert at least the first two characters of the patient record number or the last name into the search field. In order to select a record, simply click on it. It will then be highlighted in red as the record you have selected.

The screenshot shows a dark-themed search interface. At the top is a search bar containing the text '12'. To the right of the search bar are a magnifying glass icon and a close (X) icon. Below the search bar, a dropdown menu is open, displaying a patient record. The record is highlighted in red. The record details are: a female symbol, 'ATBM_04_Dermoskopie, Manuela', '3/17/1993', and 'ID: 123'. A small grid icon is visible at the bottom right of the dropdown menu.

The patient record that was most recently selected remains active in the search field of the dashboard. This selection can be cleared by clicking the x-symbol in the upper right corner.


The screenshot shows a dark-themed search bar with the placeholder text 'Search for patient name or patient ID here'. To the right of the search bar are a magnifying glass icon and a close (X) icon.

3.2.2. Making patient data anonymous

A patient's name can be made anonymous in the *Patient List* so that it is not displayed while in the presence of a third party. This previously installed setting can be modified in *Settings* area of the *Dashboard*. Anonymous patient names can, nevertheless, be selected while using the query function. (cf. 8.1)

3.2.3. Editing patient information

Should you require editing an existing patient record, simply open the Patient Information tab on the right after selecting a record from the Patient List (cf. 3.2.1). You can then modify and add details as desired. Make sure to click the Save button afterwards.

Thursday, 13 June 2013


ATBM_04 Dermoskopie, Manuela 3/17/1993
ID: 123

Patients (1st 29 patients) Patient Information


Last Name	First Name	Patient Record Number	Middle Name	Date of birth
444	444	333		16.05.2013
ATBM	Demo	1415		06.11.2012
ATBM	Puppe	435		01.08.1995
ATBM_04_Bodyscan	Manuela	146		17.01.1993
ATBM_04_Dermoskopie	Manuela	123		17.01.1993

3.2.4. Deleting a patient record

If you want to delete a record, select the record from the database first (cf. 3.2.1). Then click on *Delete Patient*.

As a security measure, you will be asked if you are sure about deleting the patient.

After confirmation, the patient record and all related images will be removed irrevocably!

Once you have successfully deleted a patient record, go back to the dashboard by clicking back.


4. Section "Total Body Mapping" (TBM)

In order to capture and process images in the Total Body Mapping (TBM) section, open it from the dashboard. We use the section when working with the automated photography.


To be able to run the section, a patient record has to be selected, which you can either create prior to executing the TBM section (cf. 3.1) or select from your patients database (cf. 3.2.1). Then click on the *Total Body Mapping* icon on the dashboard.

Should you not have provided information for the patient's height yet, you will be asked to do so now. Make use of the slide control to indicate the height in centimeters. After clicking *OK*, the info will be automatically saved with the patient information in the *Patient Management (C)* section. Here, you can still amend the info later on, if necessary.

4.1. The interface

The TBM desktop consists of six main sections:

A Menu Bar

B Patient Data Field

C Timeline

D Body Map

E Preview Window with zoom control

F SlideGallery


These sections provide the functions described below.

A Menu bar.

The menu bar is located at the top of the screen.

It is made up by the following functions:

Click *Back* to switch to the previous screen. When you click *Back* on the desktop, you will leave the module and return to the dashboard.

Baseline images are the images that are captured in an initial session. They are by default set up as reference for future screenings.

Follow up images are the images you capture following an initial session. They can be compared to the baseline images in detail thanks to the functions of the Total Body Mapping section.


Use the *Marker* button to highlight the position of a lesion on a patient's skin. The markers are numbered consecutively.


Compare images from different sessions directly next to each other in detail.

Report allows you to create a document in PDF format by default, designed specifically to match Total Body Mapping.

The *Export* function allows you to export selected images in jpeg format to local or external storage devices.

For optimal detection of malignant skin lesions, refer to the *Dermoscopy* button following a Total Body Mapping directly and take capture of dermoscopy images.

Edit image offers an easy tool for measurement. Markers can also be configured.

B Patient Data Field.

On the upper right hand the Patient Data Field indicates patient name, date of birth, patient record reference, and the number of images stored for this patient.


C Timeline.

Below the menu bar a timeline is provided, which allows you to jump between different imaging sessions. The selected session is marked as a colored square. Move the mouse cursor over an individual session and the date of the session will appear.

D Body map.

On the left hand side of the screen, four body maps guide you through the image capturing process. Positions with existing baseline images are shown in color while grey body parts still need to be captured. The current image, which is also shown in the preview window, is highlighted in color.

Each of the four body maps is divided into four segments: head and shoulder (highlighted in color above), torso, thighs, and legs. The software automatically creates a whole-body profile of four consecutive body segment images for four positions (front, right, back, and left).

Additionally, there are four positions for capturing the hand from the outside, the palm of the hand as well as both of the foot soles. A complete sequence therefore consists of 20 pictures. As each of the four sides is captured in an automatic sequence, you have to press *Capture* only eight times (cf. 4.2).

The aim of Total Body Mapping is to map the entire body.

You can individually select one out of twenty *Additional segments* if you want to capture an additional specific pose of your patient. They can be selected underneath the body maps.

E Preview window.

On the right hand side of the screen is the Preview Window. It shows the live image or the picture that was captured (framed in color in the SlideGallery).

There is a number of functions available, which allows you to adapt your view as needed: Hit the full resolution button (100%), click *Fit to screen* to adapt the image to the size of the preview window, or zoom into the image with the slider bar on the right. A click with the right mouse key allows you to make use of a magnifier.

F SlideGallery.

At the bottom of the screen, you find the SlideGallery.

It displays all images of the selected session. Move the mouse cursor over an individual session in the timeline and the date and time of the session will appear.

4.2. Photography

4.2.1. Taking Baseline pictures


Thanks to the unique CamControl technology the Photo camera Canon EOS is connected to the PC at all times and controlled through the software. This allows you to take standardized high-resolution images.

All the photo camera settings like exposure, flash, etc. are optimized by the program.


Make sure that the photo camera is connected correctly. Turn on the photo camera by pressing the *ON/OFF* switch on the PolFlash face and by pressing the On/Off switch on the photo camera. If the photo camera is not connected correctly, you will receive an error message. Also verify that the photo camera menu wheel is on *M* ("manual") before switching into the imaging mode. If this is not the case you will receive an error message.

Switch on the Laser Liner and align the Positioning mat against the Laser Liner beam. Use the poster with example pictures to position the patient as shown on the poster.

Now click on *Baseline* in the Total Body Mapping menu bar.

In the following window all positions on the body map are grey. This means that no images have been saved so far.

Start the capturing process by hitting the *Start* button.

Start

The photo camera now moves to the first position according to the body height that was indicated for the patient. In the window on the right, you will now see a live preview image and the corresponding body part will be highlighted in color on the body map.

Attention: If the tower control cannot be initialized, an error message appears. If that is the case, contact the authorized representative of FotoFinder, MAXmedical, LLC at: 6 Sirenevaia Street, office 54, 141109 Shchyokovo, Moskovskaia oblast. E-mail: support@maxmedical.ru .


Please note that the live image is only meant to give an indication for finding the correct position for your exposure. By no means does it reflect the absolute image quality.

Once the patient is positioned correctly, click on *Capture*.

Capture

The program now automatically captures the images from head to toe. Subsequently, the preview window portrays a complete image that consists of the four single images.

After automatically fitting the images into one full image, the margins of separate images can be further matched by hitting the button *Manual* and shifting the breaks using your mouse.

 A rectangular button with a black border and the word "Manual" in a sans-serif font.

If the captured image does not meet your expectations, click on *Repeat* and restart the capturing process. If you are satisfied with the image taken, click on *Save*.

 A rectangular button with a black border and the word "Save" in a sans-serif font.

The software now automatically jumps to the next view. For the remaining positions, continue the capturing process according to this principle. As soon as all of the twenty standard positions are captured, a pop-up will confirm the completion of the session.

If you want to stop the automatic capturing process intermittently, click *Emergency STOP!* The photo camera will stop immediately and all images of the photo camera drive will be removed.


 A rectangular button with a black border and the text "Emergency STOP!" in a sans-serif font.

Further exposures can only be made after clicking *Back* and starting the capturing process anew.

 A rectangular button with a black border, a white left-pointing arrow, and the word "Back" in a sans-serif font below it.


4.2.2. Additional body segments

Should you need to add further images to your patient record for detailed documentation (e.g. of wounds, the armpits, the scalp or underneath the breasts), you can do so by referring to the button *Additional segments*.

 A rectangular button with a black border and the text "Additional segments" in a sans-serif font.

Please note that capturing follow-up images for additional segments requires a baseline image to be saved in the respective position as well.

To take a picture in such a specific position, select one of the 20 available saving spots for additional segments (e.g. Option1). The preview window on the right allows you to nicely position the patient and adjust the photo camera height. Simply click on the arrows (see below) to set the correct photo camera position and then capture the image.

Once an additional segment has been added to the record, the option it is stored with is highlighted in color and marked with a camera symbol (see Option 1 above).

In order to switch back to the body map, click *Body View*.


Click Back to return to the desktop.

4.2.3. Taking follow-up pictures

In order to document risk patients' skin over time and provide for the easy detection of changes, refer to the function for capturing follow-up images.

Start the capturing process by clicking on the *Follow up* icon.

A message will appear if your baseline image session does not cover all of the twenty positions. Complete the baseline session or close the message and continue with the follow-up.


All photo camera settings are automatically restored to the same settings as in the baseline image.

After clicking *Start*, on the right hand side the Ghost function places the baseline image and the live image on top of each other. This helps you to position the patient correctly. An optimal processing and analysis of the baseline and follow-up images can only be provided when the patient is positioned identically with the transparent Ghost image. The picture below shows an example of a live image and a Ghost, which yet need to be matched nicely for the exposure.


Autofocus Liveview

Hit the *Autofocus Liveview* button for a sharp live image to precisely position the patient.


In order to easily match the baseline and follow-up image, you can further adjust the transparency level of the Ghost image by using the slider bar (see arrow). The standard transparency setting is at 30% of the baseline image.

Click on *Capture* after positioning the patient correctly.

Click on *Save* if you are satisfied with the image taken.

If the image does not meet your expectations, click on *Repeat*. You can now give it another go.

After saving an image the program automatically jumps to the next view. The body maps will guide you through the capturing process. Please position the patient accordingly and continue until you have captured all positions.

In order to take follow-up images of additional segments, continue as outlined in chapter 4.2.1. All images are automatically added to the gallery.

When you are finished with the images, you can go back to the desktop by clicking *Back*.

Capturing high-quality images always depends on good preparation. Please make sure to position the patient consistently and to capture the image in front of a back- ground that is uniform and properly lit.

4.3. Image treatment and analysis

Having saved images to a patient record, you will be able to generate an individual full-body skin profile for your patient in the next step. TBM provides a number of tools to facilitate the processing and analysis of your images. These include:

- *Editing Images*

- The *Compare* function.
- the *Viewing tools* available on the TBM desktop and within the compare function.
- the *Bodyscan* function.

- the *Insert marker* function

- as well as the direct connection to the Dermoscopy section that it includes.

4.3.1. Editing images

Edit image offers various features.

Zone A: You can zoom into an image. To do so, use the slider in the upper area or use the mouse wheel. (cf. 4.3.3)

You can display the results of Bodyscan by clicking on the Bodyscan icon. (cf. 4.3.4)

New markers can be added. (cf. 4.3.5)

Calibrate activates the tools on the right side of the screen by choosing a reference value (distance) in the image to which further areas and distances refer to. (See below)

By hitting *Back* the TBM desktop will be displayed.

Zone B:

With this feature an existing marker can be marked by pulling a window over it. Activate it by pressing the left mouse button. Now they will be displayed red/active. By again scrolling over them until the cross appears you can shift or lengthen it.

Hit *Angle* to calculate the angle between chosen lesions.

Colors for all elements can be chosen. The scroll bar below determines the width of new to draw or previously marked elements.

Regarding calibration:

With this button a certain length or distance in the image can be determined or called up.

With *diameter* all naevi and lesions can be measured in size. All results refer to the chosen calibration.


Surface allows you to measure drawn areas. All results refer to the chosen calibration.


Black Bar allows you to make patients anonymous.

Attention: Black bars and revised markers will be saved permanently in the image after the following query:

4.3.2. Comparing images

The *Compare* function allows you to compare baseline images with their respective follow-up images. To start this function, click on the *Compare* button.

Select the desired shooting session in the timeline to start the comparison (see arrows).


The upper image shows a picture from an initial session (on the right) next to a picture from a follow-up setting (on the left).

4.3.3. Navigational tools

To change the body view, move the mouse to left side of the screen. Simply click on a specific segment to change the view accordingly. The current selection is marked red and the image is shown in the preview window.


The second option for changing the body view is to move the mouse to the rim of the preview window. A navigation arrow appears.


When activated by a mouse click, it displays the subsequent picture (lateral, above or below, depending on the position of mouse). The basis for orientation is the position of image within the scheme.

In the *comparison mode* move the mouse to one of the sides of the image. The white flap opens up and the following window appears which will help you to select the requested parts of the body:


Select the segment or the additional segment you are interested in by clicking. The current selection is marked red.

The navigation arrows are also available in the *image comparison*. To navigate between the images move the mouse to the requested margins and activate the arrows by clicking.

Please note: The program will always show the same position in both compared images. Changing the position in one image will always change the position in both images.

Another element which helps you to move quickly to an adjacent image is the rubber band effect. To use it zoom into the preview window or the image comparison by scrolling the mouse wheel. Hold the left mouse and move the image into the direction of the next position.


A red rubber band appears on the edge of the screen, showing how far you have moved the image. If you cross a fixed limit, the program will switch to the next image. This function is only available if the image is not displayed completely – but when you have zoomed into the image.

4.3.4. Image preview tools

Following functions are available on the TBM desktop and in the Compare Mode to view the exposures in detail:

Zoom into an image by scrolling the mouse wheel. In order to move within the image, click into the image and push the left mouse button while moving the mouse pointer at the same time.



Move the mouse across the image while pressing the right mouse button. This initiates the magnifier function and enlarges the section that is marked by the cursor.


SmartZoom enables you to zoom into images simultaneously in the Compare Mode. If the function is enabled the *Zoom Lock* icon is marked by a colored bar (see right).

The *Fit To Screen* function shows both images scaled to the optimal screen resolution. You can choose to have the image adjusted to screen *Height* or *Width*.

When *100%* is selected, images are displayed in their original resolution.

The *Marker* function allows you to select whether or not marked nevi shall be highlighted in the images (cf. 4.3.6).

Marker Ghost displays marker from previous sessions in order to copy them into the current image. (ch.4.3.6)

Move Marker enables you to change the position of a marker by keeping the left mouse button pressed. The marker must have been activated with a click previously. You can also change its direction (click on the end of it).


When *Zoom Lock* is activated in the *Compare Mode*, the functions in the menu bar are performed in both images simultaneously. In line with this, the two separate menu bars are merged when *Zoom Lock* is activated (see below).

4.3.5. Bodyscan ATBM

In addition, the program provides a handy feature for detecting new and changed lesions with the Bodyscan tool. Thanks to a complex algorithm that compares new images with images of previous sessions, the tool allows an automatic skin analysis for risk patients. Any changes that are detected are clearly highlighted.

Analysis with Bodyscan is permitted only for fully grown patients. Otherwise, correct results cannot be guaranteed. In any case, the classification of the results rests within the responsibility of the doctor since the Bodyscan is based solely on statistical analyses.

The function can be executed in the Compare Mode (4.3.5.1), when adding markers (4.3.5.2) or out of the TBM desktop (4.1).

4.3.5.1. In the Compare mode

When performed in the Compare Mode, Bodyscan allows for a direct comparison of the exposures from different sessions based on the automatic analysis it provides.

The Bodyscan symbol in the top right corner of the image (see arrows in the picture below) shows that the nevi have been analyzed. If it is displayed in color, changes or new lesions have been detected. If it is displayed in white no changes have been identified or the image is the baseline exposure.

If you click the icon once, all nevi that have been found will be circled white. Additionally, all nevi that have experienced changes since the last session will be highlighted in color.


If you click the icon twice, only those nevi will be highlighted that have changed since the last screening. The *Body Scan* icon is highlighted as the image on the right illustrates.

The Bodyscan function is by default set up in such a way that the baseline data as well as the last session are taken for reference. However, it is possible to perform an automatic comparison between any two sessions. Simply select the two sessions you would like to compare by clicking on them in the timeline (Here: session 1 and 2):


Since the background analysis automatically compares the present session with the baseline data, the *Scan now* button is available in the menubar, if you select two other sessions for analysis in the timeline.

If the shown images were NOT calculated by the current Bodyscan version the *Scan now* symbol appears green. By clicking on that button a new analysis for this image is carried out by the newer version of Bodyscan.

4.3.5.2. When adding markers

A Bodyscan can be equally performed within the window for adding markers. To do so, simply click on the *Marker* icon at the top of the desktop of the Total Body Mapping section or in the image comparison mode (cf. 4.3.6).

Above the window that is opened subsequently, there is a Bodyscan icon available (see arrow in the above image). Click the icon and activate the Bodyscan as described in the previous section to have new and changed lesions highlighted.

4.3.5.3. On the TBM desktop


You can perform a body scan in the TBM section by clicking the *Body Scan* icon above the preview window (see arrow). Hereupon, all nevi are analyzed and circled.


4.3.6. Inserting and copying a marker

In the next step you can highlight the changed lesions and suspect skin areas that you have spotted with the help of the previous tools. In order to do so, please refer to the function for inserting a marker. This tool also boasts a direct link to the image comparison mode and Dermoscopy section for further analysis and diagnosis.

You can control the display of the results by clicking the *Body Scan* icon in the menu bar. It is possible to display three different phases of automatic analysis:

When the Bodyscan function is not activated the images are displayed without any highlights. In this case the *Bodyscan* icon in the menu bar is not marked either.

The *Marker* function can be performed:

1. on the TBM desktop as well as
2. in the image comparison mode.

Hit the *Marker* button. The preview of the image is framed red in the marker mode. To highlight suspect skin areas move the cursor over the image section that you would like to mark while pressing the mouse button. As soon as you release the mouse button an arrow will be added to mark the indicated skin area. These markers are automatically numbered in consecutive order. If no desired length is defined, the system creates an arrow with minimum length.


For the different locations from which the inserting of a marker is performed, this is done as follows:

(1) On the TBM desktop

For a detailed view, zoom in on the image by scaling up the bar in the top right corner of the screen (see arrow) or scroll into it. Then add the marker.


When markers are active navigation within the enlarged image is nevertheless possible. Keep the space-bar of the keyboard pressed and relocate the mouse in the desired direction.

(2) In the image comparison mode

While there are no markers for the initial session of this patient, several markers have been added in the follow-up session.

Should you require deleting a marker, you can do so by right-clicking it and selecting *Delete Marker* in the subsequent context menu.

Markers which were created during the baseline can be transferred by using the *Marker Ghost* feature for follow-ups. This helps you to avoid marking lesions twice. To copy markers to new overview images go to *Compare*.

In the timeline on the left choose your baseline image, on the right; choose the follow-up image to which markers should be copied to. (Or vice versa)

To get a preview of all markers in the follow-up, activate the *Marker Ghost*. The markers of all previous images are now displayed in a transparent grey.

At the end of the arrow the ongoing number of marker is displayed, e.g. "58", as well as the number of session in which the image was taken, e.g. S1.

To make the transfer, click with the right mouse button on the marker's number. In the context menu click on *Copy Marker to current image*.

Now the marker is transferred with all dermoscopy images into the follow-up image.

To position the copied marker to the correct lesion or naevi, activate *Move marker* by clicking on the button. Now pull the mouse over the marker by keeping the left mouse button pressed. The marker is displayed in bright red.

Now click on the arrow line, keep the left mouse button pressed and move the marker to the desired spot. When clicking on the white spot at the end of the arrow line the direction can be adjusted by pulling in this direction.

4.3.7. Connection to Dermoscopy

For a dermoscopy examination of single lesions please refer to the Dermoscopy section. TBM provides a direct link to this section, which enables you to add dermatoscopic exposures to selected nevi and take them to further, in-depth analysis. (cf. 5.)

4.3.7.1. Adding a dermoscopy image

In order to add a dermoscopy image to a particular skin area, it has to be marked on the body mapping image. For this purpose, please refer to the *Marker* function, which was described in the previous chapter.

Subsequently, you have the option to add a dermoscopy image to the marker and switch to the Dermoscopy section. This can be done departing from different locations:


1. from the TBM desktop or
2. the image comparison mode.

Right-click the marker you would like to add a dermatoscopic exposure to. Then select *microshooting* in the subsequent context menu:

(1) On the TBM desktop

(2) In the image comparison mode

When selecting this option by clicking on it, you will directly be taken to the capturing mode of the Dermoscopy module. Please note that you will only be able to make an exposure, if your Medicam is correctly connected and if you have the licensing for the Dermoscopy module.


There is also a cross-connection to the Dermoscopy section on the TBM desktop. Simply click on the button *Dermoscopy* (right) and you will be directly taken to the desktop screen of the Dermoscopy section. (Ch. 5.)

4.3.7.2. Displaying dermoscopy images

Additionally, saved dermoscopy images can be displayed within the TBM section and referred to for a detailed analysis.

They are displayed on the TBM desktop as well as in the image comparison mode. If dermoscopy images have been added to specific locations within a patient record, the numbering of added markers (arrows) is highlighted in color (see below):


There are highlighted numberings of markers in the baseline image for the above image comparison (left). This illustrates that there are dermoscopy images saved for the respective lesions.

Such dermoscopy images can be previewed in the TBM section. You can display the previews:


1. on the TBM desktop as well as
2. in the image comparison mode.

Then click on a highlighted marker numbering, that is, a marker for which a dermoscopy image is saved. A preview for the dermoscopy images added in this location will be displayed now (see below). On the TBM desktop and in the *image comparison mode* respectively the preview looks as follows:

(1) On the TBM desktop


(2) In the image comparison mode


You can display a detailed view of the dermoscopy images by clicking on the preview. You will then be taken to the following view:

1) On the TBM desktop

2) In the image comparison mode


In order to jump back to the previous view, simply click on the zoomed dermoscopy image. Entering the detailed dermoscopy image view from the TBM desktop you also have the option to switch back to the desktop by hitting the Back to Total Body Mapping Image button.

Back To Total Body Mapping Image

4.3.7.3. Comparing dermoscopy images


The connection to Dermoscopy also provides the possibility to compare dermoscopy images in detail. To run this function, simply right-click a marker and select *micro image compare* in the subsequent context menu.

You will then be taken to the Dermoscopy section.


To jump back to the previous screen, hit the *Back* button on the top right of the screen.

5. Section "Dermoscopy"

You can start this section right on the Dashboard by clicking on this symbol:


To get more information on this section we will now give a detailed description of its essential features.

5.1. User interface


Dermoscopy's user interface is composed of 4 main sections:

- A: Menu Bar
- B: Patient and Image Data Field
- C: Preview Window
- D: SmartGallery


5.1.1. Main menu

The Menu Bar is located in the screen's upper corner.

To jump back to the previous window, hit the *Back* button.

By clicking on the arrow a drop down menu will appear. You can now choose between *Medicam 800HD* and the *Schema selection*.


The *Marker* function allows you to add arrows to the images, to indicate lesions' localizations to be able to assign these to the dermoscopy images later.


Here you can also create overview images by using Medicam 800HD.


Here you can take dermoscopy images using Medicam 800HD.

The *Compare* button allows you to compare images from a patient's different examinations.


By hitting the *Delete* button the (orange framed) image can be deleted. You have to agree to that action again before it will be finally deleted.

The *Report* button opens the printing menu.

Click on *Import* to import overview images from your hard drive or your external memory.

The *Export* button allows you to export images to internal or external memories, e.g. CD, USB flash drive, external hard drive, etc.

5.1.2. Patient data area

The patient data field is located on the left screen below the menu bar. This field consists of all relevant patient data as well as additional text fields for the *Shooting Date*, *Treatment*, *Diagnosis*, *Study*, *Localization* and *Comment*.


The *Shooting Date*, *Treatment*, and also the *Localization* will be filled in automatically while the screening takes place. All other fields are optional.

If you would like to comment on a certain image, fill in the *Comment* field which immediately appears after having taken the picture.


Patient Information	
Last Name:	ATBM_04_Bodyscan
First Name:	Manuela
Date of Birth:	17/03/1993
Image Information	
Shooting Date:	26/03/2013 14:54
Treatment:	
Diagnosis:	
Study:	
Localization:	KV
Comment:	

5.1.3. Image preview window

The Preview Window is located on the right hand side of the screen. Here you can enlarge selected overview or dermoscopy images. The numbered and colored arrows indicate that there are dermoscopy images available for a certain lesion. By clicking on the select number the dermoscopy image will appear in the gallery. (see 5.1.4. SmartGallery).

A double click on the Preview Window enlarges the dermoscopy image displayed. In the upper screen corner you can see the following Menu Bar:

By using the slide control you can enlarge an image (pure software enlargement). When using the active enlargement and clicking on the left mouse button within the image you can move the image as you like. The *Best fit* button allows you to adjust the image to the optimal screen resolution.

5.1.4. SmartGallery

The SmartGallery in the lower screen separately displays overview and dermoscopy images. If you have not taken any images for a certain patient so far no images will be displayed in the SmartGallery. For the user interface's better illustration an already documented patient is used here.

Dermoscopy images are only displayed to the respective (orange marked) overview image.

Overview images


The field *Overview image* display all of a patient's overview images. The images are sorted by the screening date. If the gallery contains a large number of images which cannot all be displayed at once, a scroll bar will appear on the right hand side. Using this scroll bar you can see all images. The active (currently selected) image will be framed in orange.

Furthermore, the selected image will be displayed in the Preview Window above. The images' size in the gallery can be adjusted by using the slide control (located on the right hand side above the gallery).


All the current image's relevant data is displayed above the gallery. This data includes the screening date and time, treatment, image type baseline session/ follow-up session as well as the resolution.

Additionally, the size of the gallery can be adjusted to meet your needs. To adjust it, scroll toward the area between the overview images and the Patient Data Field or the Preview Window. When the cursor changes, press the left mouse button, keep it pressed and pull the field in the desired direction (up or down).

Dermoscopy images

Dermoscopy images are displayed on the lower screen end.

All dermoscopy images for the current (orange framed) overview image are displayed. As soon as you choose a marker in the Preview Window this will appear in red and only the corresponding images will be displayed in the dermoscopy images gallery.


By clicking on one of the dermoscopy images in the gallery this dermoscopy image will be displayed in the Preview Window. By double clicking in the Preview Window you can see the image in high resolution. Similar to the overview images' strip here is also an orange strip located above the dermoscopy images which contains information like screening date and time, treatment, image type (baseline session/ follow-up session) and image

resolution.

Colors within the SmartGallery:

Dark purple:	TBM baseline session (also imported ones)
Pink:	TBM follow-up session (also imported ones)
Green:	dermoscopy image
Blue:	overview image created by using Medicam 800HD.


5.2. Overview pictures

The purpose of overview images is to save the nevi's accurate position whose dermoscopy images were captured.

The single nevi will be numbered consecutively to be able to assign the dermoscopy images.

5.2.1. Overview pictures with Photo camera Canon EOS

Standardized overview images of the whole body are taken with your FotoFinder bodystudio ATBM's *Total Body Mapping* Module. This is accessible from the Dashboard.

The *Total Body Mapping* and *Dermoscopy* sections together offer comprehensive precautions against skin cancer. The software instructs you to take pictures of your patient in at least 20 standardized positions.

The *Ghost* function which is integrated in *Total Body Mapping* makes it easy to take standardized follow-up images.

The captured pictures will automatically be used as Overview images for the dermoscopy examination as soon as you add the markers within the *Total Body Mapping* section. They will automatically appear framed in color in SmartGallery. (see 5.1.4).

5.2.2. Overview pictures with Medicam 800HD

In *Dermoscopy* overview images can be taken with Medicam 800HD. This is helpful if a body part with lots of nevi needs to be specially documented after the Total Body Mapping process or if a patient has only a few lesions.

If you want to take overview images with Medicam 800HD please consider the following:

- Make sure that the room is bright and that there are no large light variations. Sometimes additional light is necessary to prevent shadows on the patient's images. Do not take pictures in front of the window while this can produce overexposed images.
- For overview images we suggest using the Medicam's 800HD floodlight feature.
- Make sure to screen all relevant nevi. Therefore sometimes an additional close-up needs to be done, e.g. if there are hidden lesions or obese patients.
- Hold the video camera still while taking shots. For close-ups, the Medicam integrated spacer helps.

Take off Medicam 800HD's dermoscopy lens (D-Scope III Plus Lens) by twisting it out and select *MediCam 800HD* in the Menu Bar.

Click on *Overview* to start the imaging mode.

Now you are in the imaging mode. Please switch on the LED light of the video camera before starting a record. The button is on the camera's back panel (also see Item 4.4.4. of Technical Manual).

Now click on the video camera's shutter release button or click on the *Freeze* button. The picture will freeze afterwards.


If the image taken is blurred or if you are not satisfied you can again activate it by hitting the *Live image* button. The capture will be deleted.


By pushing the shutter release button again or when clicking on *Save* the picture will be saved.

Before saving the captured image you need to choose its localization. Localization means choosing the patient's body part which a lesion rests on. This is chosen by clicking on the mannequin's corresponding body parts.

The mannequin can be turned around by clicking on the green arrow above it.

When scrolling over the relevant positions these will be highlighted grey. By clicking on one section this will be primed red and recorded in the localization field. The section is active until another one is chosen.

You can also add further details, e.g. *Diagnosis*, *Treatment* or *Comments*.

To save the image, please, hit the shutter release button again or click on *Save*. After the capture was saved the video camera switches back to the live image and you can take the next overview image right away.


You will go back to the Desktop when clicking on the *Back* button and you can no longer take any images.

The Patient and Image Data Field now automatically displays the screening date and time, treatment and localization.

5.2.3. Overview pictures with mannequin

Instead of clinical overview images you can also use the virtual patient's corresponding detail. Select *Scheme selection* in the drop-down menu. The virtual patient will appear and using the green arrow you can change her position as desired. Click on the corresponding body part and hit on *Save* afterwards. The virtual patient's out-take will now be displayed in the gallery as overview image and can be labeled.


5.2.4. Inserting marker

In new or saved overview images you can set additional markers at any time. A marker is always necessary to take a dermoscopy image of a lesion. The marker is a link which points at the dermoscopy image (see 5.2.4)


Important: A dermoscopy image cannot be taken without prior marking. Each dermoscopy image needs to be assigned to a certain position within the overview image.

For marking please proceed as follows:

Click on the *Marker button*.

Now, while holding the left mouse button down, draw a line from the overview image to the desired lesion. As soon as you release the mouse button the line will be saved and numbered automatically.


Number the whole overview image like that. As soon as all areas are numbered switch off the image marking by clicking on the *Marker* button again.

A patient's numbering is consecutive, that is, each number is awarded only once. It is therefore impossible to have the image wrongly assigned.

Deleting markers

Micro position marker that is no longer needed can be deleted later on. To delete it, turn off the marker in the overview image by clicking on the right mouse button and selecting the *Delete* command in the following menu.

With the line color command you can select the color for the numbering in the following menu.


If you want to skip a position in all of the subsequent follow-ups, e.g. because the lesion has been removed, you can deactivate these markers. Therefore select the marker with the right mouse button and choose *Deactivate*. You can activate the marker the same way.

5.3. Dermoscopy images with Medicam 800HD

Tips regarding dermoscopy images:


- Attach a dermoscopy lens (D-Scope II Plus Lens) to the video camera and turn on the light to do the dermoscopy images.
- Before starting the screening make sure that the lens, the glass tube and/or the polarizing filter are clean. Cleaning instructions are provided in the Technical Manual (Item 2.7.).
- To take dermoscopy images without a polarizing filter always use alcohol solutions which contain 70% alcohol. The skin needs to be generously moistened. Try to prevent bubble formation and screening of hair. Shave if necessary.

- Please make sure that, if the optional polarization filter without immersion liquid is used, sufficient disinfection is put on the skin.
- To analyze a lesion using the Moleanalyzer or FotoFinder Trichoscale Pro programs, all images need to be taken with 20x enlargement. Furthermore, we suggest using immersion liquid.

Put on a D-Scope II Plus Lens on the Medicam 800HD and select *MediCam 800(HD)* in the Menu Bar.

Choose the captioned overview image for which you want to take a dermoscopy image.

Then click on *Micro*.

Now the capture mask for dermoscopy images is displayed. The upper left window shows the live image. Next to it the most recent reference image is displayed. Below the marker position is displayed in detail on the left hand side. On the right the overview image is displayed with active colored markers.

- 1 Marker position
- 2 Live image
- 3 Reference image
- 4 Overview image

If there is no reference image yet, the corresponding window will be blank.

Now turn on lights on the back panel of the video camera. Moisten the body parts that need to be screened by using alcohol containing spray and put a glass tube attached to the video camera on the lesion that needs to be screened.


If you are using the optional polarization filter no immersion liquid is required. The procedure is the same.

The *Zoom level* button allows you to select a desired enlargement. Alternatively the enlargement factor can individually be selected using the zoom button on the video camera's back.


If there are already dermoscopy images for this position the video camera automatically sets the reference image figures. Using this, an objective comparing is made possible. You can change the image's enlargement if necessary.

The video camera will use the autofocus function. If you want to analyze a nevus using the Moleanalyzer later on, the Medicam's screening needs to be enlarged twenty times. This zoom factor also needs to be used for the hair analysis using FotoFinder Trichoscale Pro.

Now hit the release or click on the *Freeze* button. The image will be stopped right away and displayed full-screen.


If the image is blurred or if you are not satisfied with it you can activate it by hitting *Live image* and thereby deleting the image.

By pushing the camera's release or by clicking on *Save* the image will be saved.

After saving the image the software selects the following number of marker. Repeat the last steps to take further images.

The *Previous* and *Next* button allow you to jump from one marker to another. After all of an overview image's markers have been called up the system automatically go to the next overview image in the chronological order.

You can go back if you want have a closer look at a recently saved image.

Or you choose *next* if a lesion has not changed since the last screening and if you consider that it unnecessary to save it once again.

When clicking *Back* you will go back to the Dashboard.

5.4. Images comparison

5.4.1. Comparing dermoscopy images

To compare the dermoscopy images you first have to call up an overview image, by clicking on the number of the desired localization.


As a result the arrow will appear red and the number orange. In the gallery the already existing dermoscopy images will be displayed in small format (micro image preview).

Click on *Compare*. Once it is activated it will be marked orange. As an alternative you first choose a dermoscopy image that is saved in the SmartGallery to enlarge it in the Preview Window and select the *Compare* function.

The following user interface will open:


The two latest images are displayed.


Above the dermoscopy images the marker's position number is displayed. Additionally, the images of every session are listed chronologically. Furthermore, information regarding localization and shooting date are displayed.

You can compare several images of one position (number) by selecting the desired image using the scroll bar located below the dermoscopy images.

The further left you scroll the older the dermoscopy images get. Finally you will see the baseline overview image. The selected dermoscopy images/markers are tagged red.


Using the "Previous" and "Next" buttons you may go to dermoscopy images of the market with the next number, without quitting the compare screen.

You may compare only the dermoscopy images that have the same number.


5.4.2. Comparing overview images

To be able to compare overview images you first of all have to call up an overview image.

Please click on *compare*.

In this section you can have a closer look on and compare all of a patient's overview images taken so far. Use the scroll bar located below the images to call up the desired overview image. Older images are further to the left, the latest one is on the far right.

Above the overview image you see the total number of overview images taken so far as well as their chronological order. Furthermore, information regarding localization and screening date are displayed.

You can go back to the User Interface of the module Dermoscopy click on *Back*.

5.4.3. Other functions in the image comparison mode

SmartZoom

This button is set as default if you select the compare section. *Zoom Lock* allows you to zoom in and switch from one to another image simultaneously so that a precise comparison is ensured. Using the button in the upper screen allows you to zoom in the shown images simultaneously.

The 1:1 button displays both images in their original resolution. Best fit will adjust both images to the screen resolution.


To move around in the zoomed in image click on it and keep the left mouse button pressed and scroll the cursor.

Without Zoom Lock both images can be zoomed in independently. Hitting this button will deactivate this function and another zoom control will appear on the screen center.

Now you can enlarge the dermoscopy images independently of each other. The Solo button will display each image individually and in full-scale. Clicking on it once again leads you to the compare function.


Live image

The compare section allows you to take a follow-up on a clinical image. Therefore select an overview image in the screen's right hand side and click on Live image.

A preview is displayed on the left hand side. By clicking on Capture and Save you can take a new image. Thanks to CamControl the video camera's settings will be adopted.

Copying

Using the Copy button you can insert compare images in other programs, e.g. Word or Power Point. Therefore an accurate screenshot is taken, i.e. using the selected zoom factor, etc.

Select the desired images and click on Copy. Alternatively, you can hit [Ctrl]+[C]. Afterwards click on the target file using the right mouse button and click on Insert or hit [Ctrl]+[V].

Above the images information regarding the screening date and time are displayed.

Moleanalyzer

The **Moleanalyzer** helps to analyze selected dermoscopy images.

5.5. Printing and reports


You are in the SmartGallery. To print marked images (which are ticked and/or displayed in the preview window) select an image in the gallery (overview image or dermoscopy image) and hit *Report*.

FotoFinder creates a PDF file immediately and asks for its storage location. You can print the PDF or electronically hand it over to your patient on a CD or a USB stick.

If you want to print one or several dermoscopy images (ticked in the gallery) the overview images which indicate the nevi will be printed automatically as well.

Illustration: Print layout with three selected dermoscopy images

If you mark one or more overview images (by ticking them in the gallery) all corresponding dermoscopy images will be printed as well.

5.6. Optional: Moleanalyzer program

The most important purpose of the Moleanalyzer program (Version: 3.0.13.146 of 28.07.2014) is to measure pigmented spots and to objectively track changes in the follow-up. Using the Moleanalyzer you can analyze single nevi or two follow-up images at the same time. By activating an overlay function changes can be demonstrated.

The system should only support the diagnosis as a "second opinion".

Open the program after you have marked the dermoscopy images which have been enlarged by a factor of 20 in your SmartGallery (marked = orange).

The program automatically starts to analyze the marked lesion. The specially developed algorithm provides reference to conspicuousness for certain parameters, e.g. surface, diameter, edging irregularity, asymmetry, tonal value, areas, etc.

A score will be calculated afterwards which helps the doctor to put a diagnosis and will be displayed in form of a bar. (white = low risk of malignancy, yellow = elevated risk, red = needs control more often) Two lesions can be analyzed on the screen at the same time. The result can be printed immediately.


Interface of the Moleanalyzer program.

You can also use the *Overlay* function and put the baseline image and the follow-up on top of each other to reveal tiny changes. To do so, click on *Overlay*.

In the following window appears:

Important: The score calculation is only possible if

- There is a melanocytic nevus,
- Image enlarged by 20 times,
- Taken with Medicam 800HD.


Therefore please make sure that the dermoscopy images are taken enlarged appropriate if you want to analyze the lesions later on. Furthermore, please make sure that not every lesion is suitable for analysis with Moleanalyzer, e.g. if it is very small nevus.

5.7. Optional: PASIvision program

PASIvision program (Version 2.0.2.51203 of 16.12.2015) enables the evaluation of the Psoriasis Area Severity Index (PASI) in patients with the pathology. By the comparison of images made with a time interval, the program enables doctor to clearly see any changes in skin condition in psoriasis patients. PASIvision program is launched from Section 4 "Dermoscopy" of the FotoFinder universe dashboard.


Interface of the PASIvision program.

← PASI analysis results.

For PASIvision, patients are photographed with the same total body photography system. The body surface area, erythema, thickness and scaling are evaluated during the analysis, and the Total Evaluated PASI Score is calculated and presented.

5.8. Optional: FotoFinder Trichoscale Pro program

Interface of the FotoFinder Trichoscale Pro program.

The FotoFinder Trichoscale Pro program (Version: 2.2.1 of 10.11.2015) is used for measuring parameters of hair growth and loss. The module automatically determines the number and density of hair in a certain area and its anagen-telogen rate. Please read the FotoFinder Trichoscale manual to know how to make preparations for the analysis.

After having all preparation steps done (cutting and dyeing the hairs), please take a dermoscopy image of the scalp with Medicam and D-Scope II Plus Lens using 20X magnification. To do so, proceed to making dermoscopy images as described in this manual in Article 5.3. "Dermoscopy images with Medicam 800HD" on Page 92.

After you have taken a dermoscopy image, go back to the desktop and select this image in the SmartGallery. You can now start the FotoFinder Trichoscale Pro program by clicking on the corresponding button.

FotoFinder Trichoscale Pro measures the following parameters: analysis area, hair number and density, anagen-telogen rates, density and percentage of vellus and grown hair, and quantity of follicular units. It also has tools for making manual corrections. The program presents analysis results via on-screen trichogram with three different views of the analysis area and advanced diagrams and pie charts.

6. Section "Screening"

Screening Mode allows you to do a time-saving skin screening on your patient by using the Medicam 800HD. Thus you can make use of its enlargement and image quality without having to save the image taken.

You can start the Screening Mode by clicking on the following button located on the Dashboard:

6.1. User interface


When the video camera is switched on you see the user interface which contains two areas, the Menu Bar and the video camera's picture on full screen.

6.2. Main menu

The Menu Bar is located on the upper screen.

By clicking the *Back* button you will go back to the previous screen.

With *Freeze* or by pushing Medicam shutter release button once the live image will freeze.

With the *Live image* button or by once again hitting the shutter release button a Live video will be activated.

With *Copy* you can copy a frozen image in the FotoFinder clipboard and add a marker later on

The Moleanalyzer helps to analyze selected dermoscopy images (20x enlarged).

By using the *enlargement factor* you can enlarge as desired. Alternatively the enlargement factor can be adjusted through the video camera's zoom button individually.

6.3. Screening procedure

- Start the Screening Mode
- Switch on the lenses light on Medicam 800HD's back and put the camera on the skin covered with antiseptic spray. We strongly recommend using alcohol containing spray.
- Pushing the release button will freeze the displayed image.
- Now you can have a closer look at lesions. The live image will again be activated after pushing the release.

If you have the add-on module Moleanalyzer you can freeze in an image which has been shot enlarged by a factor of 20 without and analyze it right away without the need to save it beforehand.

6.4. Saving images

It is not possible to save images in the Screening mode. If the screening reveals conspicuous nevi, images taken in the Screening mode can be copied to the clipboard and assigned to a marker in *Dermoscopy*.

- To do so, copy the frozen image by clicking on *Copy* in the FotoFinder clipboard.
- Exit the Screening mode, open *Dermoscopy* and open the overview image which you want to assign the dermoscopy image to.
- In the overview image click on the corresponding marker using the right mouse button and choose the *Insert Screening* command which appears in the menu.


As soon as the corresponding overview image has been selected the dermoscopy image now also appears in the Smart Gallery.

7. Section "Query"

FotoFinder universe program offers the opportunity to systematically search your database with the help of different search criteria. Thus you do not need to search through patient lists for appropriate images for a certain presentation, but you can call up the desired images fast and directly using the query function. Start this function on the Dashboard.

Have a look at the following user interfaces:

A: Menu Bar

B: Entry Mode for search criteria

C: Result preview section.

You can start different actions right from the Menu Bar. In the upper menu you can find the standard functions.


The *Back* button allows you to go back to the Dashboard.

Here you can start a new query by putting all parameters.

By clicking *Execute query* you can activate your defined new query. The previous search will be discarded if it has not been saved.

Fields for specific queries are located below the general functions. You can carry out two different queries (full-text search, extended query) and sort the results according to certain criteria (*Sorter*).

7.1. Full text search

In the *Full text search* you can enter a certain term, e.g. melanoma. Add data in the search field. Afterwards click on *Execute query*.


All images which are assigned to a patient will be displayed and can be collated afterwards (see below). All images will include the shooting date.


Tip: By scrolling over (not clicking on) the displayed date of each picture an information section will appear. This section includes the date, the patient's name and ID.

7.2. Extended search

By clicking on *Extended query* several features you can use during the query will be displayed automatically in the upper menu.

The *Extended query* offers several drop-down menus which will help you to limit the results. Patient and Image Data can be chosen.

Through the submenu you can limit your search using further parameters.

Patient parameters:

Last name
First name
Middle name
Patient record number
Date of birth
Gender
Body height
Address
City
State
Zip
Phone
Mobile
Email
Fax

Image data parameters:

Shooting date
Diagnosis
Localization
Study
Treatment
Comment
Image name
Image type id
Image type text
Image recording device
Camera profile name

Sample query:

You are searching for dermoscopy images of lesions on left hands of men younger than 50.

1. First of all select a new filter. To do so click on:

[Click here to add a new condition] .

2. Now choose a corresponding parameter using the mouse:

Patient → Gender is equal to male.

3. Another condition now narrows down the age:

Patient → Date of birth is before (precise date) → d/m/y (written dd/mm/yy).

4. Additionally, you can search for localization

Image → Localization contains → leg (written text).

Attention: At the beginning of each query you define if all three queries are conditions, at least only one or one part of the queries has to be met. The more parameters you add, the higher the hit accuracy is.

You can save your query if you would like to. Click on the Save button in the Menu Bar's upper corner. Your query will be saved as an FFQ file in the FotoFinder universe query folder.


You can open saved queries by hitting the Load button. A new window will open and you will be able to delete the files using the right mouse button.

After several queries you can repeat the last one by clicking on Last query.

If you are familiar with other database systems, e.g. SQL, you can also display this query.

7.3. Sorter function

Click on *Sorter* and on the chosen images afterwards. Afterwards drag these, with pressed left mouse button, upwards in the *Sorter*. The system will copy the images accordingly.


The buttons for the further treatment of selected images are located on the right side.

You can delete the images in the Sorter by clicking on *Delete sorter*.

Add last sorter puts all selected images from the last sorted arrangement in the current sorter.


Export sorter (jpg) all current sorted images as single JPG to the selected area on your FotoFinder server.

Through *Print sorter* you can print images for your patients.

8. Settings

Settings can be managed right on the Dashboard itself. The button to do so is located in the upper right corner.


By clicking there the following menu bar will open. There you can manage different settings. The most important ones are described in the following section.

8.1. General


Under *General* you can adjust settings regarding the patient administration, e.g. that a new ID will be awarded each time a new patient is added. To do so, please, tick off the corresponding function.

Furthermore, patient data can be shown anonymous. A patient name will be displayed using his or her initials in the *Patient List*. The search function on the Dashboard still displays the full name.


8.2. Imaging devices

Under *Image Capture Devices* you can control the different parameters of all connected photo cameras and video cameras, e.g. the Medicam white balance and the iris value.

8.3. Other devices


Under *Devices* you can determine, that the photo camera automatically moves to a parking position at eye level after leaving the program. There it must be tied. It prevents that the photo camera moves to a low position caused by its own weight and so could be damaged when moving the device. Put the check mark at the corresponding box. Whenever you leave the program there will appear a prompt which you need to confirm.


8.4. Modules

Under *Module* different options for *Dermoscopy* are given.


Adjust these to your working method, e.g. you can order overview images by localization.

Also the printing option can be changed in *Reports*. Put the check mark in *Prompt for extended options*.

Now the window *Extended options* will open as soon as you click on *Report* when executing the program.

In this dialog you can change the printer and other options, e.g. preview, output as PDF or image.

8.5. Changing user information

During the installation process you filled in your user information for the print layout. This data is now displayed in the header. You can change the company name, personal name, address and email.

Please go to:

Start → Programs → FotoFinder.Universe → Tools → SetCustomerInfo

and open the tool by double-clicking:

A window opens where you can fill in updated information. When hitting *Accept* changes will be saved and you exit the program.

FotoFinder.Universe - Customer Informations

FotoFinder

Customer information for print layouts

Company:	FotoFinder Systems
Name:	Max Mustermann
Address 1 (Street):	Industriestrasse 12
Address 2 (City):	Bad Birnbach, Germany
E-mail:	info@fotofinder.de

 Apply

When creating a PDF file, this data will now be inserted automatically.

If there are any further questions concerning the settings please contact our authorized representative, MAXmedical, LLC at: 6 Sirenevaia Street, office 54. Phone: +7 499 502 79 44, fax +7 499 502 50 36. E-mail: support@maxmedical.ru

9 Security of data

Personal data entered and stored in the database of the FotoFinder R3 program come within the purview of the regulations on the protection of personal data. When using the database at a computerized working place the user bears responsibility for the appropriate use, for the unauthorized access of any outside persons, for the integrity, preservation and reasonable annihilation of personal data. Using the videodermatoscope "FotoFinder dermoscope" in Russian Federation the user should observe the requirements of Russian regulations regarding the security of personal data.

FotoFinder.Universe - Customer Informations

FotoFinder

Customer information for print layouts

Company:	FotoFinder Systems
Name:	Max Mustermann
Address 1 (Street):	Industriestrasse 12
Address 2 (City):	Bad Birnbach, Germany
E-mail:	info@fotofinder.de

 Apply

When creating a PDF file, this data will now be inserted automatically.


If there are any further questions concerning the settings please contact our authorized representative, MAXmedical, LLC at: 6 Sirenevaia Street, office 54. Phone: +7 499 502 79 44, fax +7 499 502 50 36. E-mail: support@maxmedical.ru

9 Security of data

Personal data entered and stored in the database of the FotoFinder R3 program come within the purview of the regulations on the protection of personal data. When using the database at a computerized working place the user bears responsibility for the appropriate use, for the unauthorized access of any outside persons, for the integrity, preservation and reasonable annihilation of personal data. Using the videodermatoscope "FotoFinder dermoscope" in Russian Federation the user should observe the requirements of Russian regulations regarding the security of personal data.

Note. The manufacturer reserves the right, without prior notice, to make minor corrections in the program interface that do not change the program integrity described in this User Manual.

On behalf of the FotoFinder Systems GmbH, Germany

Julian Mayer,

Authorized Officer of FotoFinder Systems GmbH

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