
CARDIOLINE

your profession our mission

User's Manual

ar2100view

CARDIOLINE S.p.A.

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um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc - 17/09/2013

Code 36510176

All rights reserved © **Cardioline SpA.**

CARDIOLINE® is a registered trademark of **Cardioline SpA.**

No part of this book may be reproduced in any form, nor may it be stored in a retrieval system or transmitted in any form, without written permission from:

Cardioline SpA

Via De Zinis, 6 – 38011 Cavareno (TN) - Italy.

All the information included in this manual is subject to modification without notice. In case of difference between user manual and device please refer to a representative.

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tel. +39 02 94750470 fax +39 02 94750471

Contents

1. Introduction	5
Intended use	6
1.1 How to read the manual	6
1.2 Information and recommendations relating to safe use	7
1.2.1 Instructions for the Bluetooth transmission	8
1.3 The electrocardiograph	9
1.3.1 Parts, symbols and controls	11
2. Installation and initial preparation	16
2.1 Selecting the installation site	16
2.2 Option code entry	17
2.2.1 Updating the option code	17
2.3 Loading the thermal paper	18
2.4 Power supply; control and management of the rechargeable battery	18
2.4.1 Recharging the battery	18
2.5 How to switch on the electrocardiograph	19
2.6 How to switch off the electrocardiograph	19
2.6.1 Auto power off	20
3. Preparation for use: the menu	21
3.1 How to access the menu	21
3.2 Structure of the menu	22
3.3 Menu-activated operation and customization of the electrocardiograph	24
3.3.1 "Customise Mode" ☒	24
3.3.2 The "ECG archive"	25
3.3.3 "Settings"	25
3.3.4 "Use profiles"	26
3.3.5 "General"	29
3.3.6 "Password"	30
3.3.7 "Tools"	31
4. Preparing for an ECG recording	32
4.1 Connecting the patient cable	32
4.2 Preparing the patient and applying the electrodes	32
4.3 Select recording parameters: operating mode, display and print format, speed, sensitivity, filters	34
4.3.1 Operating mode ☒ ☒	34
4.4 Display and print format	35
4.4.1 Speed of recording on the display and on paper	36
4.4.2 Sensitivity of recording on the display and on paper	36
4.4.3 Recording filters	36
5. Recording of a rest ECG	37
5.1 Patient data entry ☒	37

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Tel. +39 02 94750470 - Fax +39 02 94750471

info@cardioline.it

5.2	Recording in manual mode	39
5.3	Recording in automatic mode	39
5.3.1	Automatic calculation of the ECG parameters ☒ ☒	39
5.3.2	Automatic ECG interpretation ☒ ☒	40
5.3.3	Copy of an automatic ECG recording ☒	41
5.3.4	ECG memory: saving a recording ☒	41
5.3.5	ECG memory: archive management ☒	42
5.3.6	Saving to Personal Computer archive ☒ ☒	42
5.4	Recording in ECG Autotimer mode ☒	43
5.5	Recording in "Paper Saving" mode ☒	43
5.6	Recording in "PC ECG" mode ☒ ☒	44
5.7	Defibrillation!	44
6.	Management and control of electrocardiograph functionalities	45
6.1	Disconnected electrodes, potential defibrillation	45
6.2	Batteries low or in need of recharging	45
6.3	Print system control. Out of paper	45
6.4	Status messages and error indication: description and related event	46
6.5	Troubleshooting	47
7.	Maintenance	48
7.1	Self-test	48
7.2	Replacing the thermal paper	49
7.3	How to clean the device and the electrodes	49
7.4	How to replace the batteries	49
7.5	Maintenance of the display	50
7.6	How to clean the printer	50
7.7	Periodic checks	51
7.7.1	Waste disposal	51
7.7.2	Technical information	51
8.	Technical specifications	52
8.1	Harmonized Standards Editions	53
8.2	Standard accessories supplied	55
8.3	Options	55
8.4	Accessories	55
9.	Warranty	56
10.	Electromagnetic emission	57
11.	Disposal	61

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
info@cardioline.it

1. Introduction

ar2100view is an electrocardiograph with dual power supply (mains and rechargeable internal batteries), which in the basic configuration will:

- record an ECG exam in automatic, manual and timed mode;
- real-time display of the ECG signal on the built-in graphic display ;
- reproduce the ECG signal on 210 mm paper in various formats thanks to the high resolution thermal printer: 3/ 3+1/ 3+3/ 6x2/ 12 channels;
- organize recordings thanks to the clock, date stamp and extended alphanumeric keyboard to manage user and patient data;
- set up to 4 different user profiles to tailor how the device functions with different requirements.

ar2100view can be equipped with:

- "Memory option": storage of up to 200 full ECG exams, with no need to print out immediately on paper ("paper saving" mode);
- "ECG parameters option": automatic ECG parameter measurement program;
- "ECG interpretation option": module for automatic interpretation of the ECG signal;
- "PC archive option": for saving the exam to archive stored in a personal computer running CARDIOLINE software **cubeecg**;
- "PC-ECG option": for real time display of the twelve leads on your computer screen to allow management of patient medical records and archiving of exams in digital format using **CARDIOLINE** software. The software has an optional module for automatic interpretation of the ECG signal.

The peculiarities of the above features allow the classification of the device application field within the class IIa (Rule 10 Annex IX).

For more information on available options, contact your selected dealer.

CONGRATULATIONS ON YOUR PURCHASE. Your new computerized electrocardiograph **CARDIOLINE®** has been designed and built in compliance with the applicable regulations in force at the time when **Cardioline SpA**, Cavareno (TN) - ITALY drew up this manual. **Cardioline** operates in accordance with the requirements for quality management systems defined by UNI EN ISO 13485:2004. The system is certified by Cermet Soc. Cons. a r.l. Your new electrocardiograph has also been built in compliance with the Medical Device Directive 93/42/EEC (legislative decree 24 February 1997) and is therefore marked by the relevant CE0476 mark.

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info@cardioline.it

Intended use

ar2100view is a digital electrocardiograph, able to acquire and digitalize the 12 leads simultaneously, to visualize in real time on the integrated LCD display and to print on thermal printer, to memorize 200 ECG traces (optional), to archive them and/or transmit them to the **cubeecg** system (optional) through BT receiver for wireless PC.

The device is designed for the monitoring of cardiac functions as well as for use as an diagnostic tool.

The proper use of the product is its use in a medical environment. The device should only be used by qualified and suitably trained personnel, and must comply with the instructions contained in the User Manual.

Shelf Life: 5 years.

1.1 How to read the manual

In order to ensure the **CARDIOLINE® ar2100view** is operated in a safe and correct manner, and to appreciate its ease of use and high reliability, **the user instructions must be read carefully.**

This documentation describes the functions of your electrocardiograph including those provided by all the possible "options" available. It is therefore possible that some of the functions described may not be present in the model you have purchased. For details of the options, consult the "firmware configuration" chart that accompanies each individual appliance.

 This symbol allows you to identify the functions **not** provided on all models which must be requested specifically at the time of purchase.

 This symbol allows you to identify the functional, behavioural and operational aspects that may be conditioned by the type of configuration selected during the step of "Preparation for use: the menu".

When a given key is depicted in the body of a sentence or a paragraph, press the corresponding key on the device to perform the action.

The structure of this manual allows you to approach the use of the electrocardiograph according to your level of knowledge. If you have already had experience with **CARDIOLINE®** equipment, the initial fast-track part of each paragraph will allow you to begin working immediately. In the continuation of the paragraph, on the other hand, the single aspects of operation are discussed in more depth.

The manual gives detailed information on the use of the model **ar2100view** for traditional ECG procedures, and an introduction to the use of particular functionalities involving interaction with software and a Personal Computer.

instructions on the use of the software applications for Personal Computer, consult the special online guidelines.


For further information contact directly:

CARDIOLINE® - Customer Support Centre

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e-mail: service@cardioline.it

tel. +39 02 94750484 fax: +39 02 94750472

1.2 Information and recommendations relating to safe use

- The electrical system used by the device must be in accordance with the standard in force.
- Always use the equipment according to the instructions in this manual.
- The device is equipped with a set of standard accessories. For reasons of safety, reliability and conformity with the Medical Devices Directive 93/42/EEC, use only original accessories or accessories approved by the manufacturer.
- The device is equipped with a special long-life thermal head writing system, which allows maximum writing precision. To avoid frequent and costly replacements and repairs, always use the original paper or paper approved by the manufacturer. The manufacturer will not accept liability for any damage to the device or any other adverse effect caused by the use of unsuitable paper.
- Do not subject the device to impact or excessive vibrations.
- Do not allow liquids to penetrate inside the device. If this should accidentally occur, have the device tested by an Authorised Assistance Centre to verify its functional efficiency before using it again.
- Make sure that the value of the supply voltage corresponds to that indicated on the data plate of the device.
- If you are using the device in connection with others, ensure that: all connections are made by skilled persons; all connections comply with safety regulations; all other devices connected respond likewise to regulations. Non-compliance with regulations can cause physical harm to the patient connected and to the person operating the device. Should it be difficult to obtain the necessary information for assessing the risk of the individual connections, apply directly to the manufacturers concerned or avoid making the connections.
- In the event of other equipment being connected directly or indirectly to the patient, check for the possible risks caused by the sum of the leakage currents on the body of the patient.
- The device is protected against defibrillation discharges in accordance with standard; to ensure that the signal is restored, use only original electrodes or

- electrodes responding to IEC and AAMI standards.
- If an electrosurgical scalpel is in use, the patient cable should be disconnected from the device.
- At all events, when defibrillators or high-frequency surgical devices are being used at the same time, it is essential to take the greatest care. If there is any doubt when such devices are in use, disconnect the patient from the electrocardiograph temporarily.
- The device recognises the impulses generated by a pacemaker and does not interfere with its operation, as prescribed by standards in use at the time of drafting this manual.
- Avoid exposing the equipment to extreme temperatures, excessive dust or dirt and very salty or damp environments; observe the ambient conditions described in detail in section 8.
- Periodically check the efficiency of all accessories and of the device itself. Contact the Authorised Assistance Centre whenever the device seems to be operating irregularly.
- To prolong the life of your device, have it checked periodically by an Authorised Assistance Centre.
- **Warning:** The electrocardiograph cannot be used for intracardial applications.
- **Warning:** It is therefore necessary before activating the equipment, to make sure of the connection to ground (normally secured by the power supply cable). If grounding of the main electrical service is not certain, do not connect the device and use it powered only by the rechargeable internal battery.
- **Warning:** do not use the device in the presence of anaesthetics or volatile gases!
- **Warning:** devices for medical applications must be used only by persons who by virtue of training or practical experience are able to ensure maximum safety and effectiveness in operation. Operators must in any event read this manual carefully and familiarise themselves with the instrument before using it on a patient.
- **Warning:** the information obtained through the use of interpretive programs or other automated diagnostic aids must be validated and countersigned by qualified medical staff.

The manufacturer will acknowledge liability for the safety, reliability and functional efficiency of the device only if:

- modifications and repairs are performed by the manufacturer or by an Authorised Assistance Centre;
- the a.c. mains power supply of the building responds to current regulations;
- the device is operated according to user instructions;
- any accessories in use are those approved by the manufacturer.

1.2.1 Instructions for the Bluetooth transmission


The device complies with the standards set in the R&TTE Directive on Radio and Telecommunications Terminal Equipment.

Nevertheless, in order to protect the device from other devices that do NOT comply with the R&TTE Directive, the following instructions must be followed:

with the abovementioned standards, we recommend to put the device as far as possible from other devices that use Bluetooth transmission.

1.3 The electrocardiograph

Front view

CARDIOLINE S.p.A.
 Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT 03153711209
 Tel. +39 02 94750470 - Fax +39 02 94750471
 info@cardioline.it


CARDIOLINE

your profession our mission

ar2100view

Back and bottom view

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 06183711209
Tel. +39 02 94750470 Fax +39 02 94750471
info@cardioline.it

User's Manual: um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc

Cardioline S.p.A.

Via De Zinis, 6
Bronzetti, 8

38011 Cavareno (TN) Italy
20129 Milano (MI) Italy

info@cardioline.it
tel. +39 02 94750470

www.cardioline.it
fax +39 02 94750471

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU


User's Manual

Cardio
Sales

1.3.1 Parts, symbols and controls

1. Keyboard

Function key	Messages & Symbols displayed / Associated LED	
 Power on/off		LED on: device connected to mains power; internal battery charging
		"full" symbol: battery charged "part empty" symbol: battery power less than 30% "empty" symbol: internal battery flat; the device must be connected to the mains power for recharging
 select start operating mode		indicated electrodes not connected or insufficient contact; saturation
 interrupt current operation; stop		
 select operating mode	Auto 	Automatic recording
	Manual 	Manual recording
	Custom 	Recording mode selected in configuration phase ("Custom mode")
 Select ECG display format	Automatic mode: formats available 3, 6, 6x2	In manual mode, print and display format selections are homogeneous CARDIOLINE S.p.A. Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN) Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano C.F. e P.IVA IT 08153711209 Tel. +39 02 94750470 - Fax +39 02 94750471 info@cardioline.it

	- Automatic mode: formats available 3/ 3+1/ 3+3/ 6x2/ 12 channels - Manual mode: formats available 3, 6, 12	In manual mode, print and display format selections are homogeneous
	- On Display: Off, M1, M2 e M3 - On paper:  M1, M2 o M3	Filter deactivated, 30Hz, 25Hz e 20Hz
	Speeds available 12,5 - 25 - 50 mm/s	
	Speeds available 5-10-12,5-25-50 mm/s	
	Wavelengths available 2,5 - 5 - 10 - 20	
	Wavelengths available 2,5 - 5 - 10 - 20	
		Displays patient data file
		Print the copy of the last automatic ECG. Data will be available until the device is switched off.

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
2. **Display:** to display the ECG signal, management of functioning operations, patient data, report data and settings.

In normal operating mode:

- Patient data area display Last and First Names.
- ECG signal display area Each lead is preceded by its name.
- Selected operating mode. In mode indicates scroll functions active.
- Heart rate.
- Selected speed of printing (value on left) and display (value on right).
- Selected sensitivity of printing (value on left) and display (value on right).
- Selected print format (value on left) and display format (value on right).
- Filters on/off indicator.

- In automatic mode, free memory space available ☒. In Autotimer mode

displays the step, the time of the step and total process time.

- Time indicator.

k) Battery status indicator.

l) Operating information and error message area. In Menu mode indicates the active keys and the presence of additional options.

3. Display management keyboard

ESC

- Cancel operation, return to previous menu.

- scroll menus and information. In main screen, scroll leads on display.

MENU

- Access to Menu.

4. Extended alphanumeric keyboard for patient data management.

- Special characters are associated with each key, and can be displayed by pressing the key in sequence.

- delete text.

- confirm.

5. Patient cable connector .

6. Paper compartment door.

7. "Reset" button  : used to re-establish normal operating conditions in the event of an error that cannot be managed using the keyboard.

8. "Mains line" connector".

9. Identification plate. Refer to this data in all communication with authorised customer support centres.

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Tel. +39 02 94750470 - Fax +39 02 94750471

info@cardioline.it

2. Installation and initial preparation

CARDIOLINE® ar2100view electrocardiograph. Suggestions are given for "selecting the installation site" and "recommendations for safe use in conformity with current statutory regulations" are indicated. Also introduced are the operations involved in preparing the electrocardiograph for use, such as "loading the thermal paper", "power supply"; "control and management of the rechargeable internal battery", "switching on and off", "the menu", "set-up".

2.1 Selecting the installation site

Your electrocardiograph complies with European directives on electromagnetic compatibility. The absence of emissions damaging to radio and telecommunications transmissions is therefore assured, as also is protection from interference emitted by other systems and equipment. Nevertheless, in order to protect the device from other devices that do not comply with the abovementioned standards, we recommend the following:

- avoid the use of mobile phones near the electrocardiograph;
- place the electrocardiograph as far as possible from electrical power lines and sources of static electricity. The ECG signal can be disturbed if the electrocardiograph is placed near sources of high voltage or electricity lines;
- avoid placing the electrocardiograph close to other diagnostic or therapeutic equipment (e.g. X-ray machines, ultrasound machines, electrically operated beds, etc.) that could be a source of excessive interference and ECG signal distortion;
- if it is impossible to position the electrocardiograph at a distance from other electrical equipment, switch the other equipment off while recording an ECG.

Also, to avoid the effect of ambient conditions when recording ECG:

- record in a room where the temperature is between 20 and 25 degrees Centigrade. This precaution prevents the patient from feeling cold, which could increase shivering and contribute to muscle tremor;
- record using the battery, disconnecting the device from the mains power supply. This avoids presence of mains power disturbance of the recorded ECG signal.

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03153701209
Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

2.2 Option code entry

When the device is switched on for the first time, an option code will be required.

	OPTION CODE	
---	-------------	--

OPTION CODE

	Esci		Ok	
--	------	--	----	--

Entering the 40 digits issued by Cardioline SpA, the device will enable for the user the choice of language and options. This feature is also useful to add future options.

2.2.1 Updating the option code

Press the MENU button , select by means of the arrows , confirm submenus TOOLS, AUTOTEST, SERVICE, when the access code is requested press in

sequence: FILTER -> WIDELIGHTS -> SPEED -> AUTO and confirm .

The option code feature can be found in the service section there. Select option

code by means of the arrows  and confirm .

After entering the new code all data stored in the device will be lost and new language or options can be selected.

CARDIOLINE S.p.A.
 Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT 03153711209
 Tel. +39 02 94750470 - Fax +39 02 94750471
 info@cardioline.it

2.3 Loading the thermal paper

CARDIOLINE® ar2100view can reproduce the ECG signal on thermal paper in Z-fold packs. Thermal paper in rolls can also be used. It's necessary to plan the format page as "Letter", refer to the following paragraphs for the menus of configuration. To correctly load the different types of paper please follow the instructions below.

- Open the paper compartment.
- Prepare a new pack and position it in the compartment. Check that the red mark on the paper is on the upper left of the pack.
- Position the paper, centring it between the two paper guides. Close the cover, positioning the paper between the rubber roller and the device case.

Caution: For safety reasons, the equipment must be used only with the cover of the paper compartment duly closed.

Caution: use only original thermal paper or paper approved by the manufacturer. The use of paper that does not respond to the manufacturer's specifications could jeopardise the correct operation of the device.

2.4 Power supply; control and management of the rechargeable battery

Your electrocardiograph uses a dual power supply system: a.c. mains and a rechargeable NiMH lead battery.

The rechargeable battery is housed inside a 5+5 compartment placed on the bottom of the device and is protected against short circuits.

Caution: before using the device, it is necessary to recharge the battery completely!

Before connecting the electrocardiograph to the a.c. supply with the cable supplied, check that the mains voltage is the right voltage for the device.

Caution: the internal battery is I (first) class.

Caution: when the device is connected to the mains, the batteries are recharged automatically, even during use.

To gain maximum benefit from the characteristics of the dual power supply system follow the indications given below.

2.4.1 Recharging the battery

The battery must be recharged when the power indicator symbol is partially empty

: the reserve charge is lower than 30%.

Connect the electrocardiograph to the mains: LED • ~ lit. Complete recharging of the battery requires at least 18 hours.

For longer life, the battery should be allowed to run down and recharged completely at least every two months.

A complete recharge allows the recording of up to 300 complete ECGs (automatic recording mode, 3 channel print format, speed 25 mm/s).

Note: recording with different print sizes the maximum number of examinations recordable may differ.

If the battery should be completely rundown (symbol ) , it is still possible to make an ECG recording by connecting the device to the a.c. mains supply.

Warning: remove the batteries from the unit if the device is not used for a long period.

Warning: do not dispose of a spent battery as ordinary refuse or litter. Only use original batteries supplied by the manufacturer. The average life of the batteries is about one year.

Warning: the device must be connected only to mains with earth connection, realised following the current regulations.

2.5 How to switch on the electrocardiograph

- for at least 2 seconds.
- The display lights up indicating the active "Use Profile" When fully on the ECG signal and selected settings are displayed: ready to proceed.

Warning: if the symbol  is displayed, internal power is insufficient and the battery must therefore be recharged by connecting the device to the mains (see section 2.4). The battery will recharge even if the device is in use.

2.6 How to switch off the electrocardiograph

- The display turns off. Settings related to the last recording are maintained in memory. To know the effect of the turning off on the last automatic recording refers to section 5.3.3.

Warning: switching off is not enabled in the cases specified in the section 0. In

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2.6.1 Auto power off

The procedure is activated only after the current operation has been completed.

- After 30 min. have elapsed without any key being pressed: reserve power > 30%.
- After 1 min. has elapsed without any key being pressed: reserve power between 15 and 30%.
- After 10 sec. have elapsed without any key being pressed: reserve power < 15%.

If the auto power off function is activated, last ECG recorded, if any, and the relative settings, will be held in the memory.

Auto power off is inhibited if:

- a recording is in progress in ECG autotimer mode;
- a recording is in progress in PC-ECG mode;
- during the self-test routine;
- during "set-up";
- the device is powered by mains.

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3.1 How to access the menu

MENU

- to access the menu from the main screen.

05:26	Main Menu	
-------	-----------	--

Personalize Mode


ECG Archive

Settings

Tools

	Exit	-	Menu	
--	------	---	------	--

- to scroll through the menu items.

- executes the action associated with the symbol  displayed: access the lower level menu, select and confirm.

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info@cardioline.it

3.2 Structure of the menu

The menu is organised in four sections: "Custom mode"; "ECG archive"; "Settings" and "Tools". The following tree layout of the menu illustrates the different levels of exploration possible, and the features that can be selected. Details on the single items are given in the following headings.

- *Custom mode*
 - Paper Saving
 - PC ECG ☒
 - ECG autotimer
- *ECG Archive*
 - View
 - ◆ PC Archive
 - ◆ Print
 - ◆ Delete
 - PC Archive
 - Print List
 - Clear
- *Settings*
 - Use Profiles
 - ◆ Choose Profile
 - ◆ Select profile
 - Profile1
 - Profile2
 - Profile3
 - Profile4
 - ◆ Set Profile (Change / New / Delete)
 - User data
 - Profile name
 - Institution - User
 - Ward/Dept. name
 - Print Report 3+R
 - Auto 3+R
 - Simultaneous
 - Sequential
 - FP Rhythm Leads
 - Channel 1:
 - Channel 2:
 - Channel 3:
 - Auto print format
 - ECG only
 - ECG measures
 - Standard report

CARDIOLINE S.p.A.
 Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT 031753711209
 Tel. +39 02 94750470 Fax +39 02 94750471
 info@cardioline.it

- Extended parameters
- ECG Interpretation
 - Report standard
 - Extended settings
 - Extended interpretation
- Configure Copy
 - Automatic copy
 - No interpretation
 - No extended page
- Archive Management
 - Save
 - Automatic
 - Request
 - None
 - Deleting
 - Manual
 - Automatic
 - Autotrx last ECG
 - Yes
 - No
- ECG Autotimer
 - No. Intervals :
 - Duration of intervals:
 - Print Format
 - Auto
 - 3 leads
 - channel 1:
 - channel 2:
 - channel 3:
 - 6 leads
 - channel 1:
 - channel 2:
 - channel 3:
 - channel 4:
 - channel 5:
 - channel 6:
- Print delay
 - Yes
 - No

– General

- ♦ Date \ Time
- ♦ Measure units
 - Cm, Kg
 - Inch, Pound
- ♦ Mains filter

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT 03 53711209
 Tel. +39 02 94750470 - Fax +39 02 94750471
 info@cardioline.it

- 50 Hz
 - 60 Hz
 - No
- ◆ Display
 - Brightness
- ◆ Keyboard
 - ◆ Sensitivity
- ◆ Bluetooth
 - On
 - Off
- ◆ Recorder ID
 - Recorder ID
- Password
 - New Password
- Tools
 - About
 - Print Configuration
 - Default Configuration
 - Self-test
 - ◆ User
 - Display
 - Keyboard
 - Printer
 - Memory
 - About
 - ◆ Service
 - Only available for Technical Service

3.3 Menu-activated operation and customization of the electrocardiograph

Listed below are the operating and configuration details associated with the single items of the menu. To operate the menu refer to section 3.1.

3.3.1 "Customise Mode"

The menu allows one of the enabled modes to be added to the direct selection of

traditional operating modes (Automatic and Manual). The mode selected displayed in the lower left corner of the main screen.

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
www.cardioline.it

fax +39 02 94750471

Options / Actions available	Description
☒ Paper saving 	The ECG trace is recorded and saved without any hard copy of the signal generated. The steps and the quality of the recording are monitored through messages in the display. Thereafter, the recording can be printed or saved to PC archive. Feature associated with "memory" option.
☒ PC-ECG 	The twelve ECG leads are displayed in real time on your Computer screen where it is possible to perform all the ECG recording operations thanks to the software.
Autotimer ECG 	ECGs are recorded automatically at user-defined intervals in the "Settings" menu.

3.3.2 The "ECG archive"

The menu allows the main archive management functions to be operated.

Options / Actions available	Description
View	Displays the ECG list contained in the archive. Scroll the list and select an exam/patient. The user can now: - Archive the selected exam to a Personal Computer ("ArchPC"); - Print the selected example ("Print"); - Delete the exam selected ("Delete").  To look for a patient (research key: Last Name) in the archive.
PC archive	Transfers and saves all stored exams to PC.
Print list	Prints the list of exams currently in the memory.
Clear	Deletes all exams held in the memory. A confirmation message is displayed.
Find	To look for a patient (research key: Last Name) in the archive.

3.3.3 "Settings"

The "Settings" menu allows different users to set up the functions of the electrocardiograph as best suits the individual working methods. Configurable characteristics are grouped together in three menus: "Use profiles", "General" and "Passwords".

3.3.4 "Use profiles"

Up to four different Use profile can be saved. This facility enabled the recorder to be modified for different use requirements in just a few seconds.

Accessing the "Use Profiles" menu the user can:

- recall an existing profile: "Select Profile".
- The list of saved Profiles is displayed. Select the desired Profile; the recorder settings are changed correspondingly;
- Set a new profile: "Set Profile". A new Profile ("New") can be created, modified

("Modify"), deleted ("Delete").
confirm.

to select the desired option,

OK

The parameters that can be associated to each profile are listed below.

Options / Actions available	Description
User data	A Profile Name (8 characters), User Institution Name (30 characters), Ward/Dept. Name (16 characters) can be entered. The User Institution and Ward/Dept. Name will be printed on all documents and sent to the PC when data is saved. The profile name is used to retrieve the settings.
Patient data management	Patient data can be entered and updated at any time The following data items are present in the file: ▪ Patient ID ▪ Last Name ▪ First Name ▪ Date of Birth ▪ Sex ▪ Last Name 2 ▪ Blood Pressure ▪ Height ▪ Weight ▪ Type ▪ Medications ▪ Notes ▪ Ward/Dept. Name

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Tel. +39 02 94750470 Fax +39 02 94750471

info@cardioline.it

Options / Actions available	Description
Print Report 3+R	▪ "Auto 3+R" The way the signal is represented may also be selected when a "3+1" format is selected. Two options are available: 1. "simultaneous": all 12 leads of the 3x4 part have the same time reference; 1) "sequential": the 12 leads of the 3x4 part are represented in time sequence. ▪ "Rhythm Leads" It is possible to set the 3 reference leads.
Lead Sequence	The sequence of leads printed can be selected. As well as modifying the order in which the leads are printed out, the selection also influences the rhythm lead selection menu. Three options are available: 1) Standard; 2) Cabrera; 3) Custom: a special sequence can be set.


CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via Elli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA 03153711209
 Tel. +39 02 94750470 - Fax +39 02 94750471
 info@cardioline.it

Options / Actions available	Description
☒ Auto print format	The type of processing of the ECG trace acquired in automatic mode can be selected. The choice influences the type of document printed. Two configuration menus are available, linked to the type of automatic processing available: "ECG only", "ECG measures" and "ECG interpretation".  Inside the single menus to enable/disable a section. ▪ "ECG measures" – Standard report: represents the minimum report and cannot be de-selected. The following are reported: date and time of recording, patient data, note field, main ECG parameters calculated (heart rate; rhythm type; P, QT, QTc, PQ, QRS and QTr wave amplitudes; frontal vectors; axes). – Extended parameters: printout of table of ST depression values on the twelve leads and of the templates relative to the twelve leads. – "ECG Interpretation" – Standard report: represents the minimum report and cannot be de-selected. The following are reported: date and time of recording, patient data, note field, main ECG parameters calculated and interpretation. – Extended parameters: printout of table of ST depression values on the twelve leads and of the cycles relative to the twelve leads. – Extended interpretation: printout of table of the cycles representing the twelve leads and detailed ECG interpretation.
Configure Copy	 The "Copy" key can be configured. – Automatic print: the last recording is fully reprinted automatically. – No interpretation: of the last recording only the ECG is printed. – No extended page: of the last recording only the analysis is printed (if available).

Options / Actions available	Description
Archive management	The way in which ECG traces are saved can be set, and the memory status checked. If the device is being used in "Paper Saving" mode, saving is always automatic. ▪ "Save" - "Automatic": trace saved automatically after acquisition without action by operator. The operation is indicated on the display by the following messages: "Saving...." and "ECG in archive". - "On Demand": on completion of the recording or reporting phases, "Editor" enabled, the option "Save&Exit". Press  to continue. If you do not want to save the trace, select "Confirm&Esc". A copy of the trace remains available until a new recording is made or the device is switched off. - "None" saving disabled. ▪ "Delete" - "Manual": files can be deleted using the delete function on the Look in ECG archive menu. - "Automatic": archive files are deleted automatically after successful transmission to a PC with the ArchPC function. ▪ "Auto-transmission last ECG" - "Yes": last stored ECG will be automatically sent. The ECG is stored independently from the selection of storage and it is automatically deleted once the transmission is carried out. In case of a transmission error, the ECG remains on file for subsequent manual transmission. This applies only if: storage = no. If the archive is full, a message will be displayed. - "No": no auto-transmission
ECG Autotimer	The following can be set: number of recordings (intervals), up to 99; duration of intervals (i.e. time between recordings), up to 99 minutes; the print delay (the timed prints and those activated manually between two intervals refer to a sample acquired 10 sec before); number of leads being printed: twelve (in format set), three or six (selectable).
Print Delay	It is possible to enable/disable 10 s. print delay, referred to a real time event.

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C.F. e P.IVA IT 02153711209

Tel. +39 02 94750470 Fax +39 02 94750471

info@cardioline.it www.cardioline.it

tel.+39 02 94750470 fax +39 02 94750471

29/64

Options / Actions available	Description
Date / Time	Accessing the date/time programming masks. Use number keys for settings.
Measure units	Set units to be used for management of numerical patient data (Weight and Height). Two options are available: cm, Kg; Inch, Pound.
a.c. mains filter	To ensure the device operates correctly select the mains frequency in the area of use. Three options are available: 50 Hz, 60 Hz and none, which corresponds to a 50Hz or 60Hz filter or filter deactivated.
Display	"Brightness" may be set according to the environmental conditions of use using the buttons   .
Keyboard	Levels selectable: from level 1 (the least sensitive) to 3 (the most sensitive).
Buzzer	A Buzzer Key can be selected (a sound is heard when a key is pressed on the keyboard). The beep can be activated and deactivated using the buttons   .
Recorder ID	A numerical id (5 characters) for the recorder being used can be entered.

3.3.6 "Password"

A password (maximum 8 characters) can be entered, to be requested when reports are being prepared ("Editor" function), to access the archive (if present), to set/change Use Profiles or to restore the factory configuration (default). If you do not wish to use the protection system, leave the "Password" field blank.

Warning: to change or disable a Password you need to know the current active Password. It is therefore suggested that you carefully store the protection key used. If you lose it, contact the nearest Service Centre.

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www.cardioline.it
fax +39 02 94750471

3.3.7 "Tools"

The "Tools" menu allows the user to access system related information and activate the self-test and setting functions of the device.

Options / Actions available	Description
About	Displays the recorder identification data, software version installed and language.
Print Configuration	To print the general recorder configuration and use profiles entered.
Default Configuration	To restore the ECG to factory configuration, cancelling the settings selected.
Self-test	There are two self-test menus available: "User" and "Service". For details, see section 7.

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT 031753711209
 Tel. +39 02 94750470 / Fax +39 02 94750471
 info@cardioline.it

4. Preparing for an ECG recording

This section describes the preliminary operations required when recording an at-rest electrocardiogram with the **CARDIOLINE® ar2100view** electrocardiograph. In particular, indications are given for "connecting the patient cable", "preparing the patient", and "applying the electrodes". Also illustrated are the necessary procedures for choosing the correct recording parameters, such as "speed, sensitivity and activation of filters".

4.1 Connecting the patient cable

Connect the terminal plug of the patient cable to the connector  identified with the symbol , positioned on the right side of the device.

Note: to avoid breaking the patient cable, remove it from the connector gripping it by the plug, and without tugging.

Warning: the device is protected internally against defibrillation discharges; restoration of the signal is guaranteed as long as original electrodes are used. To ensure conditions of safety are always maintained, use only original accessories..

4.2 Preparing the patient and applying the electrodes

Careful preparation of the patient and correct positioning of the electrodes are fundamental in obtaining an ECG recording of high quality.

- First, make sure the patient is comfortable and relaxed and does not feel cold. The individual should lie back on a suitably large couch with arms and hands extended along the sides of the body: this will minimise the likelihood of the ECG trace being affected by muscle tremor.
- Clean the skin thoroughly with alcohol or ether at the areas where the electrodes will be placed.
- Connect each colour-coded plug of the patient cable to the respective electrode, observing the colour-position matches indicated below:

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03183211209
Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

Colour	Symbol	Electrode position
Red	R	right arm
Yellow	L	left arm
Green	F	left leg
Black	N	right leg
White – Red	C1	V1
White – Yellow	C2	V2
White – Green	C3	V3
White – Brown	C4	V4
White – Black	C5	V5
White – violet	C6	V6

- Apply a small amount of electrocardiograph conductive gel to the area of the skin that will be in contact with the electrode, spreading it carefully and evenly (this is not necessary when using disposable electrodes with built-in gel).
- The following figure shows the standard positioning of the electrodes.

Standard positioning of the electrodes

V1: on the 4th intercostal space, right parasternal;

V2: on the 4th intercostal space, left parasternal;


V3: on the 5th rib, between V2 and V4;

V4: on the 5th intercostal space, left hemiclavicular;

V5: on the left anterior axillary, same level as V4;

V6: on the left mid-axillary at the level of V4;

Peripheral electrodes: generally a few centimetres above ankles and wrists.

Warning: make certain that the conductive parts of the electrodes are not in contact one with another or with other metallic parts. In any event, silver and silver chloride electrodes are designed and manufactured in such a way as to minimise the likelihood of accidental contact between conductive parts and external metal objects. Ensure that the device is not affected by disturbances originating from the a.c. mains power supply (see section 2).

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Tel. +39 02 94750470 - Fax +39 02 94750471

info@cardioline.it

4.3 Select recording parameters: operating mode, display and print format, speed, sensitivity, filters

4.3.1 Operating mode ☒ ☒

The available recording modes (Operating mode) depend on the active configuration of the electrocardiograph.

-  to select the desired mode; the corresponding choice is displayed on the screen. The options available are: Automatic, Manual, Custom. To change Custom mode see "Settings".
- During printing, the active mode is printed on the information line of the printout.

Options / Actions available	Description
Automatic ECG	All 12 ECG leads are recorded simultaneously (10 seconds). The signal printed out refers to the same time period, and is saved.
Manual ECG	The selected leads are recorded and printed out. The recorded signal is in real time, i.e. the trace is produced simultaneously with its acquisition.
Customised	See "Custom mode" for details.

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
info@cardioline.it

4.4 Display and print format

To select the desired "display" format:

- the corresponding selection is displayed. The options available are: 3, 6 and 6x2.

- to scroll the leads on the display.

To select the desired "print" format:

- to select the desired format; the corresponding choice is displayed on the screen.

The available formats are as follows.

Format Type	Description
3	In "Automatic" mode: the twelve ECG leads are printed out in groups of three; each page presents one group. Total no. of pages depends on the settings. Additional pages if diagnosis aids (e.g. Interpretation) are available. In "Manual" mode: three leads per page continuously.
3+1	In "Automatic" mode: the twelve ECG leads are printed out in groups 2,5 sec for each lead (3x4) plus 10 sec. for the DII rhythm lead. Additional pages if diagnosis aids (e.g. Interpretation) are available.
3+3	In "Automatic" mode: the twelve ECG leads are printed out in groups 2,5 sec for each lead (3x4) plus 10 sec. for the DII, aVF and V5 rhythm leads. Additional pages if extended diagnosis aids (e.g. Interpretation) are available.
6	In "Manual" mode: 6 continuous leads
6x2	In "Automatic" mode: the twelve ECG leads are printed in groups of 6. Additional pages if diagnosis aids (e.g. Interpretation) are available.
12	In "Automatic" mode: the twelve ECG leads are printed out. In "Manual" mode: 12 continuous leads.

The selected format will be applied to all manual, automatic and auto timed twelve lead recordings.

4.4.1 Speed of recording on the display and on paper

To select the scrolling speed of the displayed signal:

- 
 the corresponding value is displayed. The options available are: 12.5 mm/s, 25 mm/s and 50 mm/s.

To select the printing speed:

- 

 the corresponding value is displayed. The options available are: 5 mm/s, 10 mm/s, 12.5 mm/s, 25 mm/s and 50 mm/s.

During printing, the paper scrolling speed is indicated in the information line.

4.4.2 Sensitivity of recording on the display and on paper

To select the sensitivity of recording on the display:

- 
 the corresponding value is displayed. The options available are: 2.5 mm/mV, 5 mm/mV, 10 mm/mV and 20 mm/mV.

To select the sensitivity of recording on printing:

- 

 the corresponding value is displayed. The options available are: 2.5 mm/mV, 5 mm/mV, 10 mm/mV and 20 mm/mV.

During printing, the paper scrolling speed is indicated in the information line.

4.4.3 Recording filters

If necessary, it is possible to activate filters capable of improving the legibility of the signal without modifying its morphology. Activation of the filters has an effect on both the displayed and the printed signal.

- 
 to activate the filters; Off, M1, M2 or M3 is associated with the corresponding symbol on the display. The filters available have been designed to reduce the effect of both mains disturbances and muscle tremor.

The special isoelectric anti drift filter (ADF) remains permanently activated.

During printing, the activated filters are indicated in the information line.

Warning: the filters of your ECG are very effective in attenuating disturbances and do not reduce the diagnostic content of the traces. Nonetheless, it is advisable to eliminate the cause of the interference and not only the visible effect on the trace (see sections 6.5 e 2).

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5. Recording of a rest ECG

5.1 Patient data entry

- 
 to display the patient file at any time.

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info@cardioline.it

Inside the patient file:

- display the "New" option. Confirming, data stored are cancelled.

MENU

- to display the "modify" option". Confirming, data stored can be modified.

- to move around inside the file and display the options associated with "Sex" and "Type".


- to enter data.

- to clear the character highlighted by the cursor.

or

- to confirm data.

ESC

- to interrupt the data entry phase with no changes.

The data entered in the patient file are, in sequence, "Patient ID", "Last Name", "First Name", "Date of birth", "Sex", "Last Name 2", "Blood Pressure", "Height", "Weight", "Type", "Medications", "Notes", "Ward/Dept Name".

The patient data are retained in memory until they are modified or the equipment is turned off.

Once entered, the first name and last name are displayed on the main screen.

Warning: new data will be saved only when the final item entered is confirmed.

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

info@cardioline.it

5.2 Recording in manual mode

Having selected "manual" mode (see section 4.3.1):

-  to start the recording. If the signal has not yet been initialised, the message "Clamp" is displayed.
-  to change the leads displayed during recording.
-  to change the leads printed during recording, with the exception of 12 leads format.
-  to interrupt the recording (stop).

Note: during a manual recording, it is possible to change the recording characteristics: speed, sensitivity, filters.

Warning: starting a manual recording cancels the last trace recorded in automatic mode.

5.3 Recording in automatic mode

An automatic recording allows the user to run available computation and analysis programs on the traces ("ECG Parameters", "ECG Interpretation" options), to obtain a copy of the recording, to save the recording ("Memory" option), and to transfer the ECG to a Computer. Having selected "Automatic" mode (see section 4.3.1):

-  to start the recording. If the signal has not yet been initialised, the message "Wait" is displayed.
- During the recording, progress messages are displayed: 1. "Acquisition ..."; 2. "Acquisition OK". The patient can now be disconnected.
-  to interrupt printing (stop). If the signal has already been saved ("Acquisition OK") it will still be possible to print a copy of the recording.

5.3.1 Automatic calculation of the ECG parameters ☒ ☒

The program for automatic measurement of ECG parameters allows to obtain a report of the main measures calculated at the end of every automatic recording.

- **Start:** automatic at the end of every recording.
- **Stop:** automatic at the end of the print of the report. The message "Analysis OK" is displayed.

The main items of information in the report are:

- computed value of the following parameters: heart rate; rhythm type; P, QT, QTc, PQ, QRS and QTr wave amplitudes, frontal vectors, axes.
- summary table of ST values relative to all twelve leads ;
- templates of all twelve leads .

Warning: if the parameters cannot be computed, the message "Analysis nd" is displayed. This situation may be due to excessive noise affecting the ECG trace or to incorrect positioning of the electrodes.

5.3.2 Automatic ECG interpretation

The automatic ECG interpreter is a function of the analysis program that can be used to obtain an evaluation at the end of each automatic recording.

- **Start:** automatic at the end of the recording.
- **Stop:** automatic at the end of printing the report. The message "Analysis OK" is displayed.

Warning: the ECG interpretation option provides for the execution of the ECG only on patient in rest conditions and in the absence of internal or external pacemakers.

Warning: the correctness of the results provided from the ECG interpretation option depends on the exact setting of the physical data of the patient (sex, age). The user is required to perform before the execution of the analysis, using the method described in the Manual of Use.

Warning: The automatic interpretation program in case of missing patient data considers for the analysis a person of 35 years old gender male. This information is also reported on the printout of the report.

Warning: if the device cannot perform the analysis due to poor signal quality the following error message appears on the display: "Analysis n.a." This situation may be due to excessive noise affecting the ECG trace or to incorrect positioning of the electrodes.

Note: The ECG interpreter program is structured in four parts

- 1) processing and filtering of the electrocardiographic signal
- 2) identification of the waveform and positioning of the markers
- 3) calculation of the characteristic parameters of the QRST complex
- 4) processing of the diagnosis and analysis of the rhythm.

The part of the program "processing of the diagnosis and analysis of the rhythm" provides the evaluation of the trace, and specifically:

- a) identification of the parameters that deviate from standard (these parameters are identified in the final document by an asterisk), for example, duration of

wave, lengthening of PQ interval, widening of QRS. Furthermore, the data related to the rhythm is analysed and evaluated and the related indications are given, for example sinus arrhythmia, sinus rhythm with extra ventricular systole with compensatory pause, etc. All diagnostic indications described are defined category B, in accordance with American College of Cardiology conventions.

- b) analysis of repolarisation changes, as internal or external, and the degree of intensity as reflected by variations in the ST-T segment. These diagnostic indications are defined category C.
- c) suggestions diagnostic of category A. These suggestions include classifications such as eg. infarction, left ventricular hypertrophy, right ventricular hypertrophy, or, last but not least, traces below standard. In the EKG program, the category A diagnostic suggestions are obtained using a multivariate alternative classification, that is a combination of statistical analysis and a ramified structure decisional technique. By virtue of its characteristics and the reliability of the results provided, the EKG program is in conformity with the requirements of IEC standards pertinent to programs for the automatic ECG trace evaluation.

5.3.3 Copy of an automatic ECG recording

La An ECG recorded in automatic mode, and computed ECG parameters if any, are automatically saved and can therefore be reproduced on paper any number of times.

-  to start printing a copy. The message "Copying ..." is displayed.
-  to interrupt printing (stop).

If the memory does not contain valid data, the message "No data" is displayed.

Note: the parameters can be modified before printing: filters, speed, sensitivity and print format.


Warning:

Saving modes depend on setup. If the device turns off, it is necessary to access the archive in order to obtain a copy of the last recording.

5.3.4 ECG memory: saving a recording

At the end of the automatic recording the ECG may be saved in the memory. The archive can contain up to 200 recordings for subsequent processing or transfer to PC ("PC Archive option").

- **Start:** automatic, or on demand at the end of a recording (see section 3.3.3).
- If "On Demand", the patient file is displayed at the end of an automatic recording. Two options are available "Save&Esc" and "Confirm&Esc". Select "Save&Esc" to save the trace, or "Confirm&Esc" if you do not wish to save the trace. In this case the will remain available for copies until a new recording is

made or until the device is switched off.

Warning: the number of tests in the memory and the percentage of space available can be seen on the main screen. The message "Memory full" is displayed when the memory space is completely full. After finishing the current operation, free space in the memory by cancelling exams or transferring the archive to the PC.

5.3.5 ECG memory: archive management

After selecting "ECG archive":

-  to select the desired action: 1. Display; 2. Archive PC; 3. Print List; 4. Empty.
 -  to scroll pages, if more than 1.
 -  to activate the find function (research key: last name).
 -  to confirm the desired action.
- See section 3.3.2 for description of the single operations.

5.3.6 Saving to Personal Computer archive

Your electrocardiograph is capable of transferring stored ECG traces to a Personal Computer equipped with CARDIOLINE management software. For details on the use of the application software, consult the specific manual.
The user can either proceed as described above, if the intention is to transfer the complete ECG archive.

- A detailed management of the ECG transmission can be done by using the "ECG archive" menu. Select the "ECG archive" menu and proceed as appropriate. For details of the options available, see section 3.3.2.

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03153711209
Tel. +39 02 94750470 - fax +39 02 94750471
info@cardioline.it

5.4 Recording in ECG Autotimer mode

Warning: To obtain a correct recording in "ECG Autotimer" mode, ensure that the battery and paper load are enough to perform the selected acquisition.

Operating in Autotimer mode, the device can make timed recordings of the 12 standard leads or of a group of three user-selectable leads.


After having selected "ECG Autotimer" mode:


- to start the recording the message "Paper, Battery!" will appear. If the signal has not yet been initialised, the message "Wait ..." will be displayed.
- The programmed recordings are in real time. The printouts shown information, alongside the date / time field, relating to the number of the interval and the time of recording (format: #xx yy min) as well as the caption "print from memory" if such option has been selected.

- to scroll the leads on the display during the test.

- to start manual recording inside an interval, the print format corresponds to the format selected on the display.

- to pause manual printing.

- to interrupt the test.

5.5 Recording in "Paper Saving" mode

If "Paper Saving" mode is selected the device performs and stores an automatic ECG recording including any ECG calculation and analysis programs (option "ECG Parameters", "ECG Interpretation") without any printout.

After selecting "Paper Saving" mode (see section 4.3.1):

- to start the recording. If the signal has not yet been initialised, the message "Wait ..." will be displayed.
- The patient data entry procedure now starts (see section 5.1). Patient data are mandatory.
- During the recording, progress messages are displayed: 1. "Acquisition ..."; 2. "Wait...." while the analysis is being performed (the patient can now be disconnected) 3. "Save...".

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C.F. e P.IVA IT 00153711209

tel. +39 02 94750470 - fax +39 02 94750471

info@cardioline.it - www.cardioline.it

tel. +39 02 94750470 - fax +39 02 94750471

User's Manual: um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc

Cardioline S.p.A.
Sales offices:

Via De Zinis, 6
Via F.lli Bronzetti, 8

38011 Cavareno (TN) Italy
20129 Milano (MI) Italy

- When saving is complete, copies of the recording can be printed, or the recordings can be transferred to the PC, using the "Copy" function or the "ECG Archive" menu.

5.6 Recording in "PC ECG" mode

Associated with the **Cardioline "cube"** software, your **ar2100view** becomes a PC-based acquisition system. For details on the use of the application software, consult the specific manual.

Having selected "PC ECG" mode (see section 4.3.1):

-  to start transmission. The message "Wait..." -> "Connected!" is displayed. During the transmission the ECG leads are also displayed on the recorder display.
-  to end transmission. The message "End!" is displayed.

5.7 Defibrillation!

If defibrillation occurs the symbol  is displayed. Within 10 seconds of the discharge the signal is automatically restored.

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Tel. +39 02 94750470 - Fax +39 02 94750471

info@cardioline.it

6. Management and control of electrocardiograph functionalities

6.1 Disconnected electrodes, potential defibrillation

Saturation events are controlled and monitored by your electrocardiograph. The response of the electrocardiograph depends on the current phase of operation.

Stop phase

- Symbol  displayed: critical electrode contact. The user can proceed with the recording; the "critical" electrodes are indicated in the print report in the information line (for example  indicates the critical nature of the right arm).

Manual recording phase


The event is indicated as for the stop phase. When normal conditions have been restored, the signal is centred.

Automatic recording phase

If the event is detected during acquisition (10s), the ECG is stopped and automatically returns to the stop phase. If the signal is already buffered, the print continues without interruption. The event is indicated in the same way as for all other phases.

Defibrillation

The symbol  is displayed. See section 5.7.

6.2 Batteries low or in need of recharging

When the  symbol is displayed the battery must be recharged: the reserve charge is lower than 30%.

6.3 Print system control. Out of paper

Sensors verify the correct closure of the paper compartment cover and indicate when the thermal paper is depleted. During a recording, printout is inhibited automatically and the messages "Out of paper" or "Printer!" are displayed for 3 seconds approx.

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Tel. +39 02 94750470 / Fax +39 02 94750471

info@cardioline.it

6.4 Status messages and error indication: description and related event

Listed below are the various error messages displayed and/or printed on paper when abnormal events occur. Each message is correlated to a specific condition or phase of operation.

Messages

Message	Status / Event description
Not Available!	Function or action not available for the selected operating mode
Analysis na	Automatic ECG analysis cannot be performed due to excessive signal noise
No data!	Impossible to obtain copy of the last recording
	▪ Critical or disconnected electrodes or potential defibrillator. ▪ Attention
Out of paper!	Paper finished, insert a new pack / roll
Printer!	Paper compartment cover open or not properly closed
	Batteries low

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT 0315371209
 Tel. +39 02 94750470 - Fax +39 02 94750471
 info@cardioline.it

User's Manual: um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc

Cardioline S.p.A.
Sales offices:

Via De Zinis, 6
 Via F.lli Bronzetti, 8

38011 Cavareno (TN) Italy
 20129 Milano (MI) Italy

info@cardioline.it
 tel. +39 02 94750470

www.cardioline.it
 fax +39 02 94750471

6.5 Troubleshooting

The following table summarises certain problems that may occur and the relative causes.

Problem	Cause
Isoelectric line drift	- Use of electrodes other than originals Use of electrodes in saturation - Insufficient electrode/skin contact - Electrode surface dirty - Patient moving
Interference from a.c. mains supply	- Voltage generator too close; presence of other clinical instruments (e.g. X-rays, etc.) - Patient in contact with metallic parts or with other persons
Muscle tremors	- Patient not relaxed - Peripheral electrodes adhering too tightly
Irregular paper transport	- End of paper roll - Paper roll incorrectly positioned - Use of non-original paper
Analysis impossible	Signal too unstable or noisy
No copy of trace	- Wrong copy set up - Recording interrupted before 10 sec "Real time" mode: - Manual - Autotimer - With the memory option on, if the device is switched off after the last acquisition
Abnormal signal	- Defective patient cable - Defective electrodes

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 C.F. e P.IVA IT03153711209
 Tel. +39 02 94750470 - Fax +39 02 94750471
 info@cardioline.it

7. Maintenance

In order to ensure correct and long-lasting operation of the device, it is necessary to perform the maintenance procedures listed below as often as required. In any case, if a component or an accessory needs to be replaced due to wear and tear or breaking, it is necessary to take immediate action and replace it.

7.1 Self-test

Run the User self-test procedure periodically. This performs a routine check on the functional efficiency of the display, the keys, the writing system and the memory. The user can also print out identifying information relative to the individual device. In the event of error messages being displayed, contact the **CARDIOLINE** Authorised Assistance Centre, and a technician will investigate and eliminate the causes of the trouble.

The self-test menu is accessed by selecting "Tools" -> "Self-test" -> "User". Before running the self-test procedure, ensure that there is paper loaded..

-  to select the type of test required.
-  to start the test.
- Tests available:
 - Display: pixel scan. The presence of blank areas signifies faulty operation of the display.
 - Keyboard: the position of the single keys is simulated in the display. Pressing a given key, the corresponding area of the display is energised. A lack of response in any one area indicates that the relative key is faulty.
 - Printer: the writing system generates two triangular waves, the character set in the memory, and signals with different speeds and sensitivities. Irregularities of the printing system are detectable in non-continuous lines.
 - Memory: a message relating to the status of the memory is printed.
 - About. The following items of information are printed: model identification, serial number of the device, details of software, version and language code.

Warning: do not run the service self-test without a qualified technician attendance.

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info@cardioline.it

7.2 Replacing the thermal paper

When the thermal paper is depleted, the device stops and any attempt to start recording is inhibited (see section 6.3).

7.3 How to clean the device and the electrodes

To clean the device, use a cloth moistened with water or denatured ethyl alcohol. Do not use other chemical products or household detergents.

For the electrodes: remove the electrodes from the patient cable and wash under running water. Do not scratch the electrodes and do not wash the leads box and the patient socket.

Note: the device cannot be sterilised! The electrodes can be sterilised with ethylene oxide.

7.4 How to replace the batteries

The batteries must be replaced if, not using the recorder connected to the network, the battery life is insufficient for use. The average life of a battery pack is about a year.

Warning: Only use original batteries with protection!

Warning: Before you remove or insert the batteries from the recorder make sure that the device is switched off and that the power cord is disconnected from the network; non-compliance with this warning could cause the patient or the operator an electric shock.

1. Open the battery compartment.

• •

• •

CARDIOLINE S.p.A.
Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03153711209
Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

2. After paying attention to the connection polarity of the batteries (red +, black -),

- remove the batteries from the compartment.
3. Insert the new batteries and reconnect the connector with the right polarity (red +, black -).
 4. Close the battery compartment.
 5. Before using the recorder, refer to the charging instructions described in Section 2.4.1.
 6. Reset the date and time in the menu.

Warning: Failure to observe the connection polarity of the battery pack can cause serious damage to the recorder as well as voiding the warranty.

Warning: do not discharge the batteries replaced as normal litter. Only use original batteries supplied by the manufacturer.

Note: Removing the battery pack for less than three days does not result in data loss.

Note: In case of downtime for an extended period, the battery should be removed.

7.5 Maintenance of the display

- Keep the display dry, avoiding condensate formation. If condensate forms the liquid crystals could be damaged.
- Use a soft cloth to clean the surface of the display. Use of rough cloths could scratch the surface.
- Do not use solvent-based chemical products.

Warning: if the display panel should break, do not swallow the liquid released. If this should come into contact with skin or clothing, wash immediately using copious quantities of soap and water.

High pressure during use can cause fonts to be displayed abnormally. Normal operations will be re-established.

7.6 How to clean the printer

To preserve the printer in good conditions it is necessary to regularly clean the components of the printer:

- After turning off the device and disconnect it from the mains, open the paper compartment;
- Clean the thermal head with a dry cloth without putting too much pressure on it;
- Wipe the rubber surface of the paper feed roller with a cloth soaked in alcohol;
- Wipe with a dry cloth.

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www.cardioline.it
fax +39 02 94750471

7.7 Periodic checks

To ensure correct and long-lasting operation of the device, it is necessary to carry out the following checks:

Operation	Frequency in months
Fully run down and recharge the battery	2
Cleaning the print head	3
Visual inspection of the device, the patient cable and accessories	6
Run Self-test	6
Checking and cleaning the roll paper holder	12
Check paper speed	12
Check keyboard	12
Cleaning the paper compartment and the paper sensor	12
Check patient cable and electrodes	12
Replace batteries	12

A general assessment of operability and leakage current by an authorized Cardioline must be carried out every 2 years.

7.7.1 Waste disposal

In compliance with the WEEE 2002/96/EC Directive, the product and all of his parts (patient cables, batteries, etc.) cannot be disposed of as domestic waste. For further information on the disposal, contact the competent office of your city or a Cardioline Assistance Centre.

7.7.2 Technical information

Cardioline SpA undertakes, when requested by qualified persons, to furnish the list of components used in the device, and the information necessary to repair the parts of the device considered serviceable.

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03150711209
Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

8. Technical specifications

A.c. mains power supply	Internal battery 90-250V 47/63Hz
Maximum current absorbed	500 mA a 115 V ~ ±10%
	300 mA a 230 V ~ ±10%
Mains protection	Fuse: T 1,6 A 250 V
Internal power source	rechargeable NiMH battery 10x1,2 Vdc; 2000 mAh
Int. power supply protection ...	Fuse Pico SHF SLO-BLO T 5 A Littelfuse
Classification (EN60601-1)	device class I
Applied part	Type CF
Class (Regulation 93/42/EEC) .IIa	
Defibrillation protection	Internal
Input dynamic	± 300 mV @ 0 Hz. ± 10 mV in pass band
Input impedance	>100 MΩ on each electrode
Common mode rejection	>94 dB with balanced electrode impedance
Frequency response	0,05 ÷ 190 Hz (-3dB)
Time constant	3,3 s
Acquisition	12 bit; 1000 samples/s/channel printing; 500 samples/s/channel in calculation and filters;
	Sensitivity (LSB) 5 µV/bit
Leads	12 leads in Standard, Cabrera with quality check of the connection of electrodes.
Signal memory	10 seconds for each lead in automatic
Recording sensitivity	manual and automatic: 2,5 - 5 - 10 - 20 mm/mV ± 5%
Writing system	Thermal printer a 8 dot/mm
	Usable print height 210 mm
Print format	3 / 6 / 12 / D3+1 / D3+3 / D6x2 / D12
Paper transport speed	5 - 10 - 12,5 - 25 - 50 mm/s
Thermal paper	Z-Fold pack: page 210x280mm *200FF gridded
Pacemaker recognition	Recognises pulse in accordance with current IEC standards
Filters	<i>Mains interference:</i> Modified digital notch filter 50-60 Hz, linear phase
	<i>Anti-drift:</i> Digital high-pass 0.5 Hz, linear phase, always enabled. Muscle tremor : 3 cut frequency: 20, 25 e 35 Hz).
Serial interface	Bluetooth class I
Keyboard	47 keys
Display	View 12 ECG channels, functional parameters, HR (30 - 300), including electrode connection signal
	Multi-colour TFT LCD graphic 480x272 pixel, display area 95 x 54 mm, Backlit with led (4,3")

Interpretation program	Parameter calculation (optional) ECG interpretation (optional)
Type of use	Continuous
Operating modes	manual: acquisition and printing in real time automatic: simultaneous acquisition PC-ECG: real time acquisition with display at PC (optional) Paper Saving: acquisition and saving without printing
Battery capacity	Internal battery: 80 min. with 6 channel print
Recharging time	Internal battery: 10 hours 100%
Housing protection category	IP 20
Dimensions	325 x 80 x 345 mm (length x height x depth)
Weight	3700 grams without paper
Environment conditions for use	Temperature: from +10°C to +40°C Relative humidity: from 25% to 95% (non condensing) Atmospheric pressure: from 700hPa to 1060 hPa
Environment conditions for transport and storage	Temperature: from -10°C to +40°C Relative humidity: from 10% to 95% (non condensing) Atmospheric pressure: from 500 to 1060 hPa

8.1 Harmonized Standards Editions

EN 60601-1: 1990
 EN 60601-1/A1: 1993
 EN 60601-1/A2: 1995
 EN 60601-1/A13: 1996
 Medical electrical equipment - Part 1: General requirements for safety
 EN 60601-1-1: 2001
 Medical electrical equipment - Part 1-1: General requirements for safety - Collateral standard:
 Safety requirements for medical electrical systems
 EN 60601-1-2: 2001
 Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard:
 Electromagnetic compatibility - Requirements and tests
 EN 60601-1-4: 1996
 EN 60601-1-4/A1: 1999
 Medical electrical equipment - Part 1-4: General requirements for safety - Collateral standard:
 Programmable electrical medical systems
 EN 60601-2-25: 1995
 EN 60601-2-25/A1
 Medical electrical equipment - Part 2-25: Particular requirements for the safety of
 electrocardiographs
 EN 60601-2-51: 2003
 Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential

User's Manual: [um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc](#)

Cardioline S.p.A.
Sales offices:

Via De Zinis, 6
 Via F.lli Bronzetti, 8

38011 Cavareno (TN) Italy
 20129 Milano (MI) Italy

CARDIOLINE S.p.A. 53/64
 Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
 Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
 tel. +39 02 94750470 fax +39 02 94750471
 C.F. e P.IVA n. 03537910209
 Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

performance, of recording and analysing single channel and multichannel electrocardiographs

EN ISO 13485: 2003
Medical devices – Quality management systems - Requirements for regulatory purposes

EN ISO 14971: 2009 Medical devices - Application of risk management to medical devices

EN 1041: 2008
Medical electrical equipment - Information supplied by the manufacturer of medical devices.

EN ISO 14155-1: 2009
Clinical investigation of medical devices for human subjects. General requirements.

EN ISO 10993-1: 2009 Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process

ISO 7000: 2004 Graphical symbols for use on equipment

EN 980: 2008 Graphical symbols for use in the labelling of medical devices

ETSI EN 301 489-17
v 1.2.1 (2002-08)

ETSI EN 301 489-17
v 2.1.1 (2009-05)

EMC and Radio spectrum. Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering essential requirements under article 3.2 of R&TTE Direct.

ETSI EN 300 328 v
1.7.1 (2006-10)

EMC and Radio spectrum. Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering essential requirements under article 3.2 of R&TTE Direct.

NOTE:

- ☐ The list of standard relate above, include harmonized standards for directives 93/42/EEC and 2007/47/EC (Medical Devices), for directive 2004/108/CE (EMC) and harmonized standards for directive 1999/5/EC (R&TTE). This last directive has been introduced in this document because the risk analysis of a medical equipment with Bluetooth technology inside (option of the Version of the devices ar SERIES), must include compliance to R&TTE directive.
- ☐ in this document, whenever we make mention of EN60601-1, we include the amendments EN 60601-1/A1, EN 60601-1/A2 and EN 60601-1/A13
- ☐ in this document, whenever we make mention of EN60601-1-2, we include the amendment EN 60601-1-2/A1
- ☐ in this document, whenever we make mention of EN60601-1-4, we include the amendment EN 60601-1-4/A1
- ☐ in this document, whenever we make mention of EN60601-2-25, we include the amendment EN 60601-2-25/A1

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03153711209
Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

User's Manual: um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc

Cardioline S.p.A.
Sales offices:

Via De Zinis, 6
Via F.lli Bronzetti, 8

38011 Cavareno (TN) Italy
20129 Milano (MI) Italy

info@cardioline.it
tel.+39 02 94750470

www.cardioline.it
fax +39 02 94750471

8.2 Standard accessories supplied

Code	Description	Qty
63030105	Set of 4 periph. ECG electrodes clamp Ag/Agcl	1
63030163	Set of 6 chest ECG electr.suction type Ag/Agcl	1
63050025	ECG patient cable IEC, 10 lead, plug 4mm	1
65090025	Dust cover	1
66020002	Gel tube for ECG electrodes , 260 ml	1
66020028	ECG ruler	1
	ECG Z-fold paper, 210x280mm x 200 sheets	1
	User manual	1

8.3 Options

Code	Description
9RT210SE	Real-time PC ECG (without BT adapter)
9IN210SE	ECG Interpretation
9ME210SE	ECG Memory Upgrade
9PA210SE	ECG Measurement Upgrade
9PC210SE	PC Archiving Upgrade (without BT adapter)

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info@cardioline.it

9. Warranty

Cardioline Spa hereby warrants that this machine is free from defects both as regards materials and processing for a period of 24 months for devices and 3 months for spare parts from the sales date. A document issued on delivery represents the proof of the sales date and shall be submitted in case of any request for intervention under warranty.

The warranty shall consist in replacing or repairing free of charge the parts that are faulty due to construction defects or defects in the materials used. The possible replacement of the machine is subject to the irrevocable opinion of the manufacturer. It is not possible to extend the warranty following any repair action. Intervention deriving from:

- tampering and carelessness by third parties, including assistance or maintenance by unauthorized personnel;
- lack of compliance with the instructions for use, unsuitable use or use different from that for which the equipment was manufactured;
- failure of the electrical supply;
- damages caused by fire, explosion, natural disasters;
- use of non-original consumables;
- transport carried out without the due precautions;
- use of software not related to the primary function of the machine;
- other events not connected with manufacturing defects are not covered by this warranty.

Unless otherwise specified, removable parts, accessories and the parts subject to natural wear due to their standard use are excluded from the warranty. By way of example, parts such as: patient cables, batteries, connecting cables, electrodes, glass parts, storage media, ink cartridges, etc.

Cardioline Spa disclaims all responsibility for any damage that may derive, directly or indirectly, to persons or things due to lack of compliance with all the provisions given in the user's manual, particularly as regards the installation, safety, use and maintenance of the machine as well as the non-functioning of the machine itself.

If repair and /or substitution of devices or spare parts is required, take them to the nearest service centre authorized by **Cardioline Spa** or send it to **Cardioline S.p.A.** establishment; material and labour are free of charge while the risks and the cost of transport shall be borne by the user.

After 24 months from the date of the devices' purchase and 3 months from the date of the spare parts' purchase, the warranty expires and service shall be carried out charging to the user the parts replaced and labour expenses at the applicable rates. Exceptions to these warranty conditions shall be valid only if expressly approved by **Cardioline Spa**.

User's Manual: [um_ar2100view_new_cardiolinespa_04_eng1_ru1.doc](#)

Cardioline S.p.A.
Sales offices:

Via De Zinis, 6
Via F.lli Bronzetti, 8

38011 Cavareno (TN) Italy
20129 Milano (MI) Italy

info@cardioline.it
tel. +39 02 94750470

www.cardioline.it
fax +39 02 94750470

10. Electromagnetic emission

Electromagnetic emission Directives and Manufacturer's declaration

- This device requires special precautions regarding electromagnetic compatibility and must be installed and used in accordance with the EMC information contained in the following document.
- Portable and mobile radio communication units (mobile telephones, pagers, etc...) can interfere with the operation of the medical device.
- Incorrect results can be determined below the minimum amplitude value of the patient's physiological signal of 0.1mV.

Compatibility table

This device must be used in the electromagnetic environment as specified below. The customer or the user must guarantee that it is used in this environment.

Emission test	Conformity	Electromagnetic environment	Emission
RF Emissions (CISPR 11)	Group 1	The device uses RF energy only for its internal operation. For this reason, its RF emissions are very low and probably will not cause any interference in nearby electronic units	
RF Emissions	Class B	The device is suitable for use in any premises, including domestic areas and where directly connected to a public, low voltage, electricity supply that supplies buildings used for domestic purposes	
Harmonic emissions EN 61000-3-2	Class A*		
Voltage fluctuation emissions/flicker EN 61000-3-3	not applicable		

* only applicable for ECG **elan** and **elan1100** models

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Tel. +39 02 94750470 - Fax +39 02 94750471

info@cardioline.it

The device is designed to operate in the electromagnetic environment specified below. The customer or the user must make sure that it is used in this environment.

Immunity test	Level indicated by EN 60601-1-2	Conformity level	Guide to the electromagnetic environment	Immunity
ESD (EN 61000-4-2)	± 6kV contact; ±8kV in air	conforms	The floor must be of wood, concrete or ceramic. If covered with synthetic material, the relative humidity must be at least 30%	
Fast Transient Burst (EN61000-4-4)	± 2kV for power supply and ; ±1kV for the signal conductors	conforms	The power supply network quality must be typical of a commercial or hospital environment	
Over voltage (EN 61000-4-5)	± 1kV differential mode; ± 2kV common mode	conforms	The power supply network quality must be typical of a commercial or hospital environment	
Voltage dips, brief voltage Interruptions and variations on the power supply line	<5% Ut (>95% dip in Uv) For 0,5 cycles 40% Ut (60% dip in Uv) For 5 cycles 70% Ut (30% dip in Uv) For 25 cycles <5% Ut (>95% dip in Uv) For 5 sec	not applicable	The device operates with buffer batteries and is designed for continuous operation	
EN 61000-4-11				
Network frequency magnetic field (50/60Hz) (EN 61000-4-8)	3 A/m	conforms	Network frequency magnetic field should have levels typical of installations in commercial or hospital environments	

CARDIOLINE S.p.A.

Sede legale: Via De Zinis, 6 - 38011 Cavareno (TN)
Sede commerciale: Via F.lli Bronzetti, 8 - 20129 Milano
C.F. e P.IVA IT 03153711009
Tel. +39 02 94750470 - Fax +39 02 94750471
info@cardioline.it

User's Manual: um_ar2100view_new_cardiolinepa_04_eng1_ru1.doc

Cardioline S.p.A.
Sales offices:

Via De Zinis, 6
Via F.lli Bronzetti, 8

38011 Cavareno (TN) Italy
20129 Milano (MI) Italy

info@cardioline.it
tel.+39 02 94750470

www.cardioline.it
fax +39 02 94750471

ed
nt.

This device is designed to operate in the electromagnetic environment specified below. The customer or the user must make sure that it is used in this environment.

y

Immunity test	Level indicated by EN 60601-1-2	Conformity level	Electromagnetic environment	Immunity
Conducted RF EN 61000-4-6	3 Veff From 150 KHz to 80 MHz	V1= 3 V rms	Portable and mobile RF communication units must not be used nearer to any part of the device, including its cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency Recommended separation distance $d = [3,5/V1] \cdot \sqrt{P}$ $d = [3,5/E1] \cdot \sqrt{P}$ $d = [24/E1] \cdot \sqrt{P}$ Where P is the max. nominal output of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended distance in metres (m). The field intensity of fixed RF transmitters, as determined by an electromagnetic investigation of the site, should be lower than the conformity level in each frequency interval. Interference can be detected near units marked with the following symbol: 	
Irradiated RF EN 61000-4-3	3 V/m From 80 to 2,5 GHz	E1= 3 V/m		

Example of separation distance with respect to transmitter power:

Max. nominal output (W)	Separation distance at transmitter frequency (m)		
	From 150 KHz to 80 Hz	From 80 to 800MHz	From 800 MHz to 2.5 GHz
0.01	0.24	0.24	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	24	24	23

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11. Disposal

Information for users:

Pursuant to art. 13 of the Legislative Decree 25th July 2005, no. 151 "Implementation of European Directives 2002/95/CE, 2002/96/CE and 2003/108/CE, concerning the disposal of waste electric and electronic devices".

The "barred waste bin" symbol printed on this medical device indicates that the product must be disposed of separately from other waste at the end of its useful life. The user must, therefore, contact the distributor or manufacturer when the equipment is to be disposed of at the end of its life. Adequate separate collection to subsequently send the equipment that is no longer in use for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and on health, and promotes the recycling of the materials of which the device is made.

The illegal disposal of the product by the user is punishable by the fines specified in the applicable national regulations.

The equipment subject to disposal have to be sent to Cardioline OOO (Russia) which, in turn, sends it to the manufacturer, namely, Cardioline S.p.A (Italy) for its disposal in the territory of EU according to the relevant directives.

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