
Medtronic

Medtronic

Medtronic Powered Surgical Solutions (MPSS)
4620 N. Beach Street
Fort Worth, Texas 76137
USA
www.medtronic.ru

tel 817.788.6400
fax 817.788.6225

Custodian:

Jenna Groves
Senior Regulatory Affairs Specialist
Medtronic Powered Surgical Solutions

28-July-2017
Date


IFUs: Midas Rex Legend Pneumatic High-Speed System; Midas Rex MR7
Pneumatic High-Speed System; Triton Power Surgical Instrument
System and Important Information for Triton Surgical Instruments
Surgical Pneumatic Instruments with Accessories

ATTESTATION


STATE OF TEXAS
TARRANT COUNTY

On this 28 day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, IFUs, issued by Medtronic Powered Surgical Solutions, presented to me by the document's custodian, Jenna Groves, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Notary:

Official Notary Signature

Subscribed and sworn to before me
This 28 day of July, 2017.

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GLOSSARY OF UNIVERSAL SYMBOLS

MIDAS REX
Legend

The following symbols may appear within this manual, on product labeling, or on the product itself:

	Attention, see Instructions for Use.	TUBE	Tube Control
R_x Only	United States federal law restricts this device to sale by or on the order of a physician.		Air Pressure Relief
REF	Reference Number		Use with
LOT	Lot Number	MULTI-USE DISPOSABLE ATTACHMENT	Multi-Use Disposable Attachment
SN	Serial Number	Instrument Case	Instrument Case
	Quantity	Lubricant/Diffuser	Lubricant/Diffuser
	For single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.	Dissecting Tool	Dissecting Tool
≈	Approximately equal to	Attachment	Attachment
STERILE	Non-Sterile	Control Unit	Control Unit
STERILE R	Sterilized by Gamma Irradiation	Refurbished	Refurbished
	Use by date	Accessory	Accessory
	Date of manufacture	REGULATOR	Regulator
	Temperature Limitations	Bone Mill	Bone Mill
	Unlock	MOTOR	Motor
	Lock	Brush	Brush
	On	Adapter	Adapter
○	Off	! USA	USA Only
	Finger-Operated Control		Manufacturer
	Foot-Operated Control		Do not dispose to unsorted municipal waste.
	Compliant with European Council Directive MDD 93/42/EEC.	TOOL	Tool Control
EC REP	Authorized Representative in the European Community		Consult Instructions for Use

Midas Rex Legend High Speed Pneumatic System

This symbol is used to alert the reader to important safety and precautionary information. When displayed on the actual device or packaging, it refers the user to the instruction manual.

The words Warning, Caution and Note carry special meanings and should be carefully reviewed.

WARNING

The personal safety of the patient, operator and/or operating room staff may be involved. Disregarding this information could result in injury to the patient, operator and/or operating room staff.

CAUTION

These instructions point out special service procedures or precautions that must be followed to avoid damaging the instrument. A damaged instrument could result in injury to the patient, operator, and/or operating room staff.

NOTE

This provides special information to make maintenance easier or important operating instructions more clear.

Read and understand warnings before use of the LEGEND SYSTEM.

Important Safety Instructions

WARNING

General

DO NOT use any other components except for Medtronic Midas Rex Legend Attachments, Legend Dissecting Tools and Legend accessories in combination with a Legend Motor. Use of devices other than Legend Systems or associated parts or components may cause injury to the patient, operator and/or operating room staff. Use of devices other than Legend Systems or associated parts or components will damage the motor and attachments and will void manufacturer's limited warranty.

DO NOT interchange components of the Legend System with components from the Medtronic Midas Rex Classic Motor System, the Midas Rex III Motor System or the Mednext® Bone Dissecting System and their associated attachments and dissecting tools.

DO NOT operate Legend System in a Magnetic Resonance Imaging Suite.

DO NOT use a Legend Motor and Legend Attachments before proper cleaning and sterilization protocol is followed.

DO NOT dissect bone without adequate irrigation to prevent thermal necrosis.

DO NOT operate the Legend Motor without eyewear protection.

DO NOT use Legend System components if damage is apparent or if components do not run properly. The Legend System must be inspected for damage prior to each use.

- Conduct a visual inspection of the motor's exhaust hose for cracks or tears.
- Check attachments for proper appearance. Install attachment and dissecting tool, then briefly run motor.
- Check motor for overheating and leaking lubricant at motor collet.
- Check attachment for overheating.
- Check dissecting tool for flail.

DO NOT attempt to use an overheating motor or attachment.

For metal transection, observe the following safety precautions:

- a. Eyewear protection is essential.
- b. Irrigate copiously to cool the cutting surfaces.
- c. Protect the wound site from metal debris.
- d. Use a clamp or grasping device to control loose fragments during transection of any metal component.

Legend Motors or Legend Attachments which fail due to extended use may as a result of such failing allow a component to detach and fall from the motor or attachment and may cause patient injury.

Legend Motors

Heavy side loads may cause overheating of the motor to the point where the motor is uncomfortable to hold.

- Never place an overheated motor on the patient or patient draping during surgery.
- Mitigate the overheating by discontinuing use and rest the motor by using intermittently, or wrap the motor/attachment interface with a moist sterile towel.
- If the motor is passed off, the receiver must grasp the motor by the proximal end close to the motor hose.

DO NOT attempt to use a motor that is leaking lubricant. Legend Motors which fail due to extended use and/or lubricant seal failure may allow foreign matter to migrate or emit from the tool collet end of the motor. Non-sterile fluid may leak into the surgical site. As a result, measures may be required to protect the patient from infection, per the physician's discretion. If a Legend Motor requires servicing or refurbishing, return motor to MPSS.

DO NOT attempt to change a dissecting tool or attachment while the motor is running, or when the motor or attachment is in an overheated state.

DO NOT use a Legend motor if the motor continues to run after releasing the foot pedal or finger control. If a Legend Gold Touch Motor or Legend Pneumatic Control Unit requires servicing or refurbishing, return to MPSS to be refurbished.

The Legend Motor's exhaust hose may have an oily film on the external surface from pressure and/or temperature differentials following sterilization. Wipe exhaust hose with a sterile cloth prior to use. After cleaning, if motor continues to have oil on the exhaust hose, return the motor to MPSS to refurbish.

DO NOT place motor, attachment and dissecting tool on the patient or in an unsecured location during surgery.

DO NOT use Legend Pneumatic Control Unit finger/foot control in the "☛" position to operate a Legend Gold or Legend Platinum Motor. The control is used in the "☛" position only for the Legend Gold Touch Motor.

DO NOT install or remove an attachment and dissecting tool from the Legend Gold Touch unless the safety slide is in the "☛" position.

DO NOT remove Legend Gold Touch Motor from Pneumatic Control Unit unless compressed gas has been released from the hoses.

Legend Attachments

DO NOT use a dissecting tool without the appropriate attachment as injury may occur to patient and/or operator.

DO NOT use excessive pressure, such as bending or prying, on attachments or dissecting tools. This may cause tool to bend or break and cause injury to patient, operator and/or operating room staff.

DO NOT use a Legend Attachment if any part of the attachment appears to be bent, loose, missing or damaged. If a Legend Attachment requires servicing or refurbishing, return attachment to MPSS.

DO NOT force or modify any attachment or dissecting tool to fit a Legend Motor.

Heavy side loads and/or long operating periods may cause overheating of the attachment to the point where the attachment is uncomfortable to hold.

- Never place an overheated attachment on the patient or patient draping during surgery.
- Mitigate the overheating by discontinuing use and rest the attachment by using intermittently, or wrap the attachment interface with a moist sterile towel.
- If the attachment is passed off, the receiver of the attachment must handle cautiously.

Remove Legend Footed Attachments cautiously and slowly per instructions to avoid injury to the operator.

Legend Dissecting Tools

DO NOT use a dissecting tool without the appropriate attachment as injury may occur to patient, operator and/or operating room staff.

DO NOT use or re-process a dissecting tool that has already been used. Using a dissecting tool that has already been used may cause injury to the patient, operator and/or operating room staff. Using a dissecting tool that has already been used may affect motor performance and will void manufacturer's limited warranty.

DO NOT use excessive pressure, such as bending or prying, on attachments or dissecting tools. This may cause tool to bend or break and cause injury to patient, operator and/or operating room staff.

DO NOT use an attachment and dissecting tool combination if flail or vibration occurs.

DO NOT use metal cutting dissecting tools on bone.

DO NOT use a twist drill or contra-angle dissecting tool at an operating pressure over 80 psi (5.5 bar).

DO NOT use a dissecting tool if the sterile packaging is damaged.

If a dissecting tool package is opened, but the tool is not used or contaminated, the tool can be re-sterilized by steam sterilization. Remove tool from all original packaging and place into an approved autoclave package. Steam sterilize as follows:

High Vacuum Steam	270° F (132° C) for 5 minutes
Gravity Displacement	270° F (132° C) for 15 minutes

The re-sterilized tool must be used promptly following re-sterilization. If rust or corrosion is encountered after re-sterilization, do not use the re-sterilized tool.

Dissecting tool flutes are sharp and can perforate surgical gloves. Always handle with care. Do not spin, twist, or pull a tool by grasping the flutes. Gently pull on dissecting tool shaft to ensure that it is fully seated and properly installed in motor.

Legend Pneumatic Control Unit

DO NOT use the Legend Pneumatic Control Unit to operate any other systems except for Legend Motors and TRITON® Power Surgical Instrument Systems.

DO NOT disassemble equipment before gas pressure is released from Legend Pneumatic Control Unit and associated hoses.

DO NOT carry Legend Pneumatic Control Unit by the N2 DISS pressure hose.

Legend Lubricant/Diffuser

Failure to properly secure Legend Lubricant/Diffuser may cause injury to operator and/or operating room staff.


During a lengthy surgical procedure with extended motor use, lubricant may drip from lubricant/diffuser.

LEGEND SYSTEM COMPONENTS

LEGEND MOTORS

LEGEND GOLD


High speed, high torque pneumatic motor used to dissect bone and biomaterials.

LEGEND GOLD TOUCH


High speed, high torque pneumatic motor with a variable speed finger control. This motor may be used to dissect bone and biomaterials.

 Instruction manual will use “” symbol to indicate information on Legend Gold Touch.

LEGEND PLATINUM

High speed, compact pneumatic motor, with a reduced sound output and slightly decreased power when compared to the Legend Gold and Gold Touch. This motor may be used to dissect bone and biomaterials.


NOTE

Each Legend Motor has a Lubricant/Diffuser housing at the end of the motor hose.

Legend motors and associated hoses do not contain latex.

 WARNING

DO NOT use any other components except for Medtronic Midas Rex Legend Attachments, Legend Dissecting Tools and Legend accessories in combination with a Legend Motor. Use of devices other than Legend Systems or associated parts or components may cause injury to the patient, operator and/or operating room staff. Use of devices other than Legend Systems or associated parts or components will damage the motor and attachments and will void manufacturer's limited warranty.

NON-DISPOSABLE LEGEND ACCESSORIES

LEGEND PNEUMATIC CONTROL UNIT

Controls speed of Legend Gold and Platinum Motors with a variable speed foot pedal, provides finger/foot control for operation of the Gold Touch Motor, supplies port for Legend motors and port for Triton Power Surgical Instrument System. N₂ DISS pressure hose is 20' long. A longer Medtronic Midas Rex N₂ DISS pressure hose is available and may be changed on-site

NOTE

Legend Pneumatic Control Unit is provided with a certified Supply Pressure gauge.

Legend Pneumatic Control Unit has a handle for ease of re-positioning unit while on floor or for carrying unit. A hole is provided in the base for hanging unit on a wall hook. Port plugs are attached to prevent debris ingress when a port is not in use.

CAUTION

- Use Legend Pneumatic Control Unit to operate only the Legend Motors and/or TRITON Power Surgical Instrument Systems.
- Legend Pneumatic Control Unit and pressure hose should not be sterilized.





LEGEND REGULATOR

Controls the delivery pressure of compressed gas to the Legend Pneumatic Control Unit. The pressure gauges monitor cylinder pressure (right gauge) and delivery pressure (left gauge).


NOTE: Outlet pressure gauge accurate to +/- 12 psi.

LEGEND INSTRUMENT CASES


Used to organize, customize and contain Legend equipment. Available in several sizes with individual labeled tray inserts.


N₂ DISS to Male Schrader Adapter

N₂ DISS to Air DISS Adapter

N₂ DISS to WF4 Adapter

Legend Motor Wrench

**NON-DISPOSABLE LEGEND
ACCESSORIES (continued)**

**LEGEND ADAPTERS
N2 DISS TO MALE SCHRADER
ADAPTER**

Allows for the Pneumatic Control Unit's N2 DISS pressure hose connection to be attached to a Female Schrader in-house gas connection.

A N2 DISS to Female Schrader Adapter is also available for connection of the Pneumatic Control Unit's N2 DISS pressure hose to a Male Schrader in-house gas connection.

N2 DISS TO AIR DISS ADAPTER

Allows for the Pneumatic Control Unit's N2 DISS pressure hose connection to be attached to an Air DISS in-house gas connection.

N2 DISS TO WF4 ADAPTER

Allows for the Pneumatic Control Unit's N2 DISS pressure hose connection to be attached to a Midas Rex Safety Valve Regulator previously used for Midas Rex Classic or Midas Rex III Motors. The in-line oiler must be removed from the Safety Valve Regulator.

LEGEND MOTOR WRENCH

Used to align arrows on motor collet flats prior to installation of a Legend attachment.

NON-DISPOSABLE LEGEND ACCESSORIES (continued)

LEGEND ATTACHMENTS

Provides support and stability to the rotating dissecting tool.

Legend motor attachments are available in various designs to facilitate a variety of surgical procedures. Attachments vary in length, diameter, and overall design. They are marked and color-coded to correspond with their associated dissecting tools. A few of the Legend attachments available are listed in the table below.

Attachment	Example	Other Details/Options
Standard Straight Attachments	AS09	
Standard Angled Attachments	AA14	
Straight Variable Exposure Attachments	AVS07	
Angled Variable Exposure Attachments	AVA07	
Fixed Footed Attachments	AF01	
Rotating Footed Attachments	AF01R	
Telescoping Attachments	AT10 (base) TT12A (tube)	The telescoping attachment requires the use of the AT10 attachment base, as well as a telescoping tube. Tubes are available in straight, curved, or hooded form.
Contra-Angle Attachment	AC16	
Metal Cutting Attachment	ASMC	
Perforator Attachment	AD01	Available in 800 RPM or 1000 RPM form.
5/32" Jacobs Chuck Attachment	AD02	
Bone Mill Attachment	BM100	

WARNING

DO NOT use any other components except for Medtronic Midas Rex Legend Attachments, Legend Dissecting Tools and Legend accessories in combination with a Legend Motor. Use of devices other than Legend Systems or associated parts or components may cause injury to the patient, operator and/or operating room staff. Use of devices other than Legend Systems or associated parts or components will damage the motor and attachments and will void manufacturer's limited warranty.

NOTE

Angled and straight attachments with the same length, marking and color band share the same dissecting tool. Curved and straight telescoping tubes with the same length, marking and color band share the same dissecting tool. Example: The 14-AM straight and 14-AM angled attachments are 14 cm long, marked 14-AM and have a green color band. All dissecting tools with the prefix 14 (14MH30) may be used in either the 14-AM straight or 14-AM angled attachment.

Color band on attachment is for visual reference. Match the nomenclature on the attachment with the nomenclature on the dissecting tool packaging.


DISPOSABLE LEGEND ACCESSORIES

LEGEND LUBRICANT/DIFFUSER

Provides lubrication to Legend Motors, diffuses spent air from motor and filters oil from exhausted air.


⚠ CAUTION

- **DO NOT** exceed 1 hour drill time per surgery for a Legend Lubricant/Diffuser.
- **DO NOT** re-use a used Legend Lubricant/Diffuser.
② Legend Lubricant/Diffuser is a single use, disposable product.
- **DO NOT** steam or EO sterilize the Legend Lubricant/Diffuser. The lubricant/diffuser will melt into the lubricant/diffuser housing with potential damage to the motor. The lubricant/diffuser must be removed from the motor's lubricant/diffuser housing prior to sterilization.


TELESCOPING TUBES

Provides support to the rotating dissecting tool. Telescoping tubes are disposable following multiple uses and should be discarded when heat or excessive vibration is noticed or insertion of tools becomes difficult.

CLEANING BRUSHES

Used to clean debris from lumen of attachments and telescoping tubes. Sized for an internal bore diameter of 3.2 mm, 2.4 mm or 1.2 mm in Legend Attachments and Telescoping Tubes.

DISPOSABLE LEGEND ACCESSORIES (*continued*)

LEGEND DISSECTING TOOLS

Legend Dissecting Tools' intended use is for cutting bone and biomaterials.

Legend Dissecting Tools are provided as sterile.

NOTE

Dissecting tool's size and geometry may contribute to excessive vibration at certain speeds. Increase or decrease speed by changing pressure setting or foot/finger control. Change to a different dissecting tool.

Dissecting tools with "L" identification are longer tools intended for light bone dissection. The increased dissecting tool stem/head configuration may affect dissection stability.

Surgeons should familiarize themselves with the performance of the dissecting tools before use.

CAUTION

- **DO NOT** use any other dissecting tool other than Legend Dissecting Tools with a Legend Motor.
- Match the nomenclature and color code on the Legend Dissecting Tool packaging to the same nomenclature and color code on the Legend Attachment.
- Use only Legend Dissecting Tools with a Legend Motor.
- ② Legend Dissecting tools are single use devices.

WARNING

DO NOT use or re-process a dissecting tool that has already been used. Using a dissecting tool that has already been used may cause injury to the patient, operator and/or operating room staff. Using a dissecting tool that has already been used may affect motor performance and will void manufacturer's limited warranty.

DO NOT use an attachment and dissecting tool combination if tool flail or excessive vibration occurs.

DO NOT use a twist drill or contra-angle dissecting tool at an operating pressure over 80 psi (5.5 bar).

NOTE

Match the nomenclature and color code on the Legend Dissecting Tool packaging to the same nomenclature and color code on the Legend Attachment.

LEGEND DISSECTING TOOL NOMENCLATURE:

9 cm attachment length for use in 9-M attachment 9MH30 3.0 mm dissecting tool head diameter ≈ M-8
Match Head dissecting tool tip design Classic legacy nomenclature (reference only)

7 cm attachment length for use in 7-6ST straight or angled attachment 7BA40-MN Mednext tip design ≈ 24B6ST
Ball dissecting tool tip design 4.0 mm dissecting tool head diameter Mednext legacy nomenclature (reference only)

8 cm attachment length for use in F2-B1 footed attachment or 8-B attachment F2/8TA23 2.3 mm dissecting tool head diameter ≈ B-1
Tapered dissecting tool tip design Classic legacy nomenclature (reference only)

8 cm attachment length for use in 8-B attachment 8TD114 4.0 mm depth stop 1.1 mm dissecting tool head diameter
Twist Drill dissecting tool tip design

DISPOSABLE LEGEND ACCESSORIES (*continued*)

ADDITIONAL DISSECTING TOOL IDENTIFICATION FOLLOWING DISSECTING TOOL HEAD DIAMETER

DIAMOND COATED DISSECTING TOOLS

Available in:

- D Diamond, e.g., 15MH30D
- DC Diamond, coarse, e.g., 10BA50DC
- DX Diamond, extra coarse, e.g., 15CY65DX

BALL DISSECTING TOOLS

Available in:

- L Long ball — elongated tool stem, e.g., 7BA20L-MN
- F Fine ball — fine fluted head, e.g., 7BA40F-MN
- C Carbide ball, e.g., 7BA40C-MN

TAPERED DISSECTING TOOLS


Available in:

- S Spiral cutting flute, e.g., F2/8TA23S

CYLINDER DISSECTING TOOLS

Available in:

- L Long cylinder — elongated cylinder head design, e.g., 14CY50L
- B Ball tip cylinder — rounded tip design on cylinder shape, e.g., 14CY50B

General Guidelines For Attachment and Tool Applications

These are general guidelines for dissecting tool applications and are not an all-inclusive listing.




SURGERY

COMMONLY USED LEGEND ATTACHMENTS



COMMONLY USED TYPES OF LEGEND DISSECTING TOOLS

SUGGESTED LEGEND MOTOR(S)

SPINE

8-B 9-M 14-AM straight or angled 15-A straight or angled	Match Head  Elongated spherical design allows controlled, delicate dissection. For laminotomy, foramenotomy, entry hole, nerve decompression, osteophyte removal...	Legend Gold Legend Gold Touch
	Ball  Helical cutting flutes dissect bone from a variety of approach angles. For debridement, decortication, corpectomy, osteophyte removal...	
	Oval  Helical cutting flutes and curved design blend acorn and ball styles to vary dissection efficiency with approach angle. For decortication, laminotomy, entry hole, nerve decompression, osteophyte removal...	
	Hole Maker/Saw  Matched sets of Hole Makers and Hole Saws are efficient and effective for interbody fusion...	
	Cylinder  Effective bone sculpturing and planing. For graft shaping, debridement corpectomy, decortication, interbody fusion, fusion takedown...	
	Acorn  Curved designs varies dissection efficiency with varied approach angles. For entry hole, laminotomy, bone shaping, debridement, corpectomy, decortication, fusion takedown...	
Telescoping	Match Head  Elongated spherical design allows controlled, delicate dissection. For laminotomy, foramenotomy, entry hole, nerve decompression, osteophyte removal...	
Footed	Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, laminectomy...	

NEUROSURGICAL – CRANIAL

7-6ST straight or angled 8-B 9-M 10-9ST straight or angled 14-AM straight or angled 15-A straight or angled	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, mastoidectomy, acoustic neuroma, skull base...	Legend Gold Legend Platinum
	Ball  Helical cutting flutes dissect bone from a variety of approach angles. For debridement, decortication, entry hole, nerve decompression, mastoidectomy, acoustic neuroma...	
	Twist Drill  Helical design with stop produces a hole with a precise depth. Ideal for plating...	
	Acorn  Curved design varies dissection efficiency with varied approach angles. For entry hole, bone shaping, debridement...	
Telescoping	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection, mastoidectomy, acoustic neuroma, skull base...	
Footed	Tapered  Slender design for precise dissection with minimal bone loss. For craniotomy, transection, osteotomy...	

General Guidelines For Attachment and Tool Applications (continued)

These are general guidelines for dissecting tool applications and are not an all-inclusive listing.

SURGERY	COMMONLY USED LEGEND ATTACHMENTS	COMMONLY USED TYPES OF LEGEND DISSECTING TOOLS	SUGGESTED LEGEND MOTOR(S)
---------	-------------------------------------	--	------------------------------

GENERAL SURGERY AND PLASTIC SURGERY
(CRANIOFACIAL/MAXILLOFACIAL/STERNOTOMY)

7-6ST straight or angled 8-B 9-M 10-9ST straight or angled 14-AM straight or angled	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection...	Legend Gold Legend Platinum
	Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, sinus dissection...	
	Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, entry hole, suture hole, midface advancement...	
	Twist Drill  Helical design with stop produces a hole with a precise depth. Ideal for plating...	

EAR, NOSE AND THROAT

7-6ST straight or angled 10-9ST straight or angled	Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, mastoidectomy, facial recess, cochlear implant bed, ossiculoplasty, labyrinthectomy, middle fossa approach, acoustic neuroma...	Legend Platinum
---	--	-----------------

ORTHOPAEDICS

8-B 9-M 14-AM straight or angled 21-TU 26-R Footed Telescoping	Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, cement removal, acetabular revision...	Legend Gold Legend Gold Touch
	Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, suture hole...	
	Acorn  Curved designs varies dissection efficiency with varied approach angles. For bone shaping, debridement, decortication, cement removal, acetabular revision...	
Footed	Cylinder  Effective bone sculpturing and planing. For graft shaping, debridement, decortication, cement removal, acetabular revision...	
	Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping...	

BIOMETALS/BIOCERAMICS/BIOMATERIALS

MC	Metal Cutter  Cutting flutes or diamond wheel design remove metals, ceramics and other biomaterials effectively from a variety of approach angles. For cutting rods, pins, plates, implants, screws...	Legend Gold Legend Gold Touch Legend Platinum
----	---	---

POWER SOURCE REQUIREMENTS

Motors	Required Operating (Dynamic) Pressure	Nominal Operating (Dynamic) Pressure	Approximate Flow Rate Required	Gas Type
Legend Gold	80 – 120 psi	100 psi	12 cubic feet/min	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar	6.9 bar	340 liters/min	
Legend Gold Touch	80 – 120 psi	100 psi	12 cubic feet/min	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar	6.9 bar	340 liters/min	
Legend Platinum	80 – 120 psi	100 psi	12 cubic feet/min	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar	6.9 bar	340 liters/min	


CAUTION

DO NOT run a Legend Motor at an operating pressure below or above the required operating pressure range. Operating pressure below 80 psi (5.5 bar) may not provide proper lubrication to the motor. Operating pressure above 120 psi (8.3 bar) may damage or reduce the life of the motor.

Pneumatic Connections

1. Regulator Hose (N2 DISS)
2. N2 DISS to Air DISS Adapter Gas Source
3. N2 DISS to Male Schrader Adapter Gas Source
4. Gas Source (N2 DISS)
5. Regulator
6. DISS/WF4 Adapter
7. Regulator

CAUTION

If you are using the Midas Rex Safety Valve Regulator instead of the Legend Regulator, you must replace the in-line oiler with the DISS/WF4 adapter before use.

1. Set Legend System non-running (static) pressure between 80–120 psi (5.5–8.3 bar) at compressed gas source. Operating (dynamic) pressure may be adjusted as needed when the motor has an attachment and dissecting tool installed and is running. (See Pre-operative Procedure #9)
2.
 - a. Remove cover from Lubricant/Diffuser.
 - b. Do not remove inner foil seal of Lubricant/Diffuser.

3. Insert Lubricant/Diffuser into lubricant/diffuser housing located at base of motor hose.

4. Rotate Lubricant/Diffuser from “” position to “” position. It will come to a complete stop with a tactile click when it is fully seated and locked into the housing.

WARNING

Failure to properly secure Legend Lubricant/Diffuser may cause injury to operator and/or operating room staff.

NOTE

The lubrication rate is pre-set by the factory and, therefore, should not be altered.

Cartridges should not be left under pressure for more than 4 hours. If this occurs, an audible air leak may be heard from the diffuser housing. If the cartridge has been pressurized for longer than 4 hours, it should be replaced.

CAUTION

- **DO NOT** use a Legend Motor without a lubricant/diffuser installed in the motor's lubricant/diffuser housing.
- **DO NOT** use a Legend Lubricant/Diffuser if it appears to be damaged or if inner foil seal is punctured.
- **DO NOT** exceed 1 hour drill time per surgery for a Legend Lubricant/Diffuser.
- **DO NOT** attempt to remove Legend Lubricant/Diffuser while the motor is in operation.
- **DO NOT** re-use a used Legend Lubricant/Diffuser.
 - ② Legend Lubricant/Diffuser is a single use, disposable product.
- **DO NOT** refill a used Legend Lubricant/Diffuser.
- **DO NOT** steam or EO sterilize the Legend Lubricant/Diffuser. The lubricant/diffuser will melt into the lubricant/diffuser housing with potential damage to the motor. The lubricant/diffuser must be removed from the motor's lubricant/diffuser housing prior to sterilization.

DO NOT sterilize. Lubricant/diffuser will melt.

5. Connect the Legend Motor hose to the Legend port on the Legend Pneumatic Control Unit.

WARNING

DO NOT pinch, kink, obstruct or step on the motor hose. This may cause the motor hose to burst with potential injury to the patient, operator and / or operating room staff.

👉 If using the Legend Gold Touch, slide the finger/foot control on the Legend Pneumatic Control Unit to the “” position. This will automatically depress and lock the foot pedal. The control will not lock into the “” position unless the motor is connected into the Legend port.

If using the TRITON Power Surgical Instrumentation System in conjunction with the Legend Motor, connect the TRITON hose into the TRITON port on the Legend Pneumatic Control Unit.


WARNING

DO NOT use Legend Pneumatic Control Unit finger/foot control in the “” position to operate a Legend Gold or Legend Platinum Motor. The control is used in the “” position only for the Legend Gold Touch Motor.


The Legend Motor's exhaust hose may have an oily film on the external surface from pressure and/or temperature differentials following sterilization. Wipe exhaust hose with a sterile cloth prior to use. After cleaning, if motor continues to have oil on the exhaust hose, return the motor to MPSS for refurbish.

6. Prior to installation of a Legend Attachment and Dissecting Tool onto a Legend Motor, ensure that arrows on the motor collet flats are in proper alignment. If the arrows are not aligned, use the Legend motor wrench to turn the rotational collet (closest to the motor case) until its arrow is aligned with the arrow on the stationary collet (farthest away from the motor case.)


NOTE

A Legend Attachment will not seat onto a Legend Motor if the arrows on the collet flats are not in alignment.

DO NOT use any other components except for the Legend Motor Wrench to align arrows on motor collet

- 👉** Before installation of attachment and dissecting tool on the Legend Gold Touch Motor, ensure that the safety slide is in the “” position.

Install a Legend Attachment and Dissecting Tool

Straight attachment installation is listed below. For more detailed information on all the various attachment and dissecting tool installation procedures, refer to the Attachment and Dissecting Tool Section.

7. Install a Legend Attachment and Dissecting Tool:

- Slide straight attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when attachment is fully seated.

- Insert dissecting tool into attachment with a slight rotational motion. A tactile and audible click is observed indicating that the tool is fully seated.

CAUTION

Match the nomenclature and color code on the Legend Dissecting Tool packaging to the same nomenclature and color code on the Legend Attachment.

- Rotate attachment to the "🔒" position on the motor case. Gently pull on dissecting tool shaft to ensure proper installation.

CAUTION

Legend Motors should only be operated when the attachment is in the "🔒" position.


- To activate the Legend Gold or Platinum Motors:
Depress the foot pedal on the Pneumatic Control Unit.
Allow motor to run for approximately 30 seconds. An increase in speed may be noticed as lubrication is provided to the motor.

WARNING

DO NOT use the Legend System if motor is leaking lubricant at motor collet/attachment.

DO NOT use the Legend System if dissecting tool is flailing.

DO NOT pinch, kink, obstruct or step on the motor hose. This may cause the motor hose to burst with potential injury to the patient, operator and / or operating room staff.

👉 To activate the Legend Gold Touch Motor:

- For the Legend Gold Touch, be certain the finger/foot control on the Legend Pneumatic Control Unit is locked in the “” position.
- With the Safety Slide in the “” position, install and lock an attachment and dissecting tool on the motor.
- Lift the Finger Control and pull the Safety Slide to the “” position.
- Depressing the Finger Control will now activate the motor. Allow motor to run for approximately 30 seconds. An increase in speed may be noticed as lubrication is provided to the motor.

WARNING

DO NOT use the Legend System if motor is leaking lubricant at motor collet/attachment.

DO NOT use the Legend System if dissecting tool is flailing.

See Troubleshooting section.

NOTE

When using longer attachments with the Legend Gold Touch, the extension lever on the Finger Control may be extended to provide better digital control.

CAUTION

When using angled attachments with the Legend Gold Touch, extend the lever on the Finger Control and hold by the attachment.

CAUTION

- If the motor has an uncomfortable temperature to touch/hold following sterilization, allow an adequate “cool down” period.
 - If the motor has an uncomfortable temperature to touch/hold and an adequate “cool down” period was used, check that a new Legend Lubricant/Diffuser is properly attached to the motor’s lubricant/diffuser housing to help cool down motor. Run the motor briefly to assist in its cool down.
9. **Operating (dynamic) pressure may be checked diagnostically at the Legend Pneumatic Control Unit while the motor is running.** Operating pressure will decrease slightly from the non-running (static) pressure setting when the motor is activated. Adjust operating pressure as needed at the compressed gas source until supply pressure gauge on Legend Pneumatic Control Unit reads within a range of 80 – 120 psi (5.5 – 8.3 bar) as required.


To ensure optimal motor performance, set nominal operating pressure to 100 psi (6.9 bar).

If the surgeon requests additional power during a procedure, compressed gas source may be turned up to an operating pressure of 120 psi (8.3 bar).


To decrease pressure, turn down the in-house compressed gas source or loosen the pressure handle on the regulator. Push down on the pressure relief at the pneumatic control unit to exhaust excess pressure in the hoses. Then re-adjust pressure as needed.

Legend Gold, Legend Gold Touch and Legend Platinum Lock

Legend Gold, Gold Touch and Platinum Motors should only be operated when the attachment is in the “” position on the motor.


Motor will operate.

Motor will not operate.

 Legend Gold Touch Lock

The Legend Gold Touch will operate with the safety slide in the “|” position and will not operate with the safety slide in the “○” position.

Motor will operate.

Motor will not operate.

In addition, the Legend Gold Touch will not operate if the finger/foot control on the pneumatic control is in the “” position, unless the foot pedal is activated also.

Legend Gold Touch Motor will not operate.

DISASSEMBLY OF EQUIPMENT

1. Depressurize System.

Turn off the tank or wall outlet compressed gas source.

Release gas from the pressure hose by depressing the pressure relief on Legend Pneumatic Control Unit.

⚠ WARNING


DO NOT disassemble equipment before gas is released from Legend Pneumatic Control Unit and associated hoses.

DO NOT remove Legend Gold Touch Motor from Pneumatic Control Unit unless compressed gas has been released from the hoses.

2. Release Legend Motor Hose from Pneumatic Control Unit by holding the motor's exhaust hose firmly and depressing Legend port connection.

When releasing the Legend Gold Touch, depressing the Legend port connection will automatically release the "☝" position on the finger/foot control to the "☞" position. The foot pedal will move upward to its non-running position.

If using the **TRITON Power Surgical Instrumentation System**, hold the TRITON blue hose firmly and remove by turning TRITON quick connect coupler clockwise.

3. Unscrew Lubricant/Diffuser from housing. Discard Lubricant/Diffuser.

NOTE

During removal of lubricant/diffuser, oil residue may be present. An absorbent cloth should be used to facilitate removal.


⚠ CAUTION

DO NOT steam or EO sterilize the Legend Lubricant/Diffuser. The lubricant/diffuser will melt into the lubricant/diffuser housing with potential damage to the motor. The lubricant/diffuser must be removed from the motor's lubricant/diffuser housing prior to sterilization.

4. Remove attachment and discard used dissecting tools in an appropriate container. Straight attachment removal is listed below. For more detailed information on the various attachment and dissecting tool removal procedures, refer to the Attachment and Dissecting Tool Section.

- To remove a straight attachment, hold motor in palm of hand. Rotate attachment to the "☞" position.
- Remove dissecting tool from the attachment and discard dissecting tool.
- Use thumb and index finger to lift attachment off motor.

⚠ WARNING

DO NOT use or re-process a dissecting tool that has already been used. Using a dissecting tool that has already been used may cause injury to the patient, operator and/or operating room staff. Using a dissecting tool that has already been used may affect motor performance and will void manufacturer's limited warranty.

STERILIZATION GUIDELINES

Legend Motor / Attachments / Tubes

Warnings and Precautions	Do not soak/submerge Legend devices. Do not use ultrasound to clean Legend devices. Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. The use of a washer-disinfector for cleaning may cause a premature degradation in performance. Allow an adequate cooling period after steam sterilization. Use ONLY nylon cleaning brushes. Non-nylon cleaning brushes leave residue that may prevent the tool from being secured properly in the handpiece.
Limitations	Verify functionality prior to re-use.

INSTRUCTIONS	
Point of Use	No particular requirements. Follow hospital procedures.
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.
Preparation for decontamination	No particular requirements. Follow hospital procedures.
Cleaning: Automated (Do NOT use ultrasonic washer)	Review the washer-disinfector warning above, before using this cleaning method. Manually rinse attachments/tubes under tap water, until no visible soil is noticed, before placing them into the automatic washer. Remove devices from instrument trays before placing into washer baskets. Orient devices following recommendations of the washer/disinfector manufacturers. After cleaning, lubricate attachments as in step 12 below. Recommended Washer Cycle Pre-Wash: Cold tap water, 2 min. Wash: 66oC, 5 min. using a neutral enzymatic detergent, pH 6.0-8.0 Rinse: Hot tap water, 1 min.
Cleaning: Manual	- Wipe all external surfaces of the motor and hose, and wipe inner surface of oiler housing with a cloth dampened with a neutral enzymatic detergent, pH 6.0-8.0. - Brush motor case and collet with a nylon brush dampened with a neutral enzymatic detergent. If using the Legend Touch motor, be sure to brush under the finger control lever. - Rinse motor thoroughly under running water, collet end pointed down. Dry collet and motor with lint free towel. - Verify that devices are visually clean after manual cleaning. - Wipe all attachments and telescoping tubes with a cloth, dampened with a surgical instrument cleaning solution. - Immerse the head of Contra-Angle attachments in surgical instrument cleaning solution and run the motor for 1 minute. - Other attachments and tubes may be mechanically agitated in cleaning solution, but not soaked or immersed. - A nylon brush dampened with a surgical instrument cleaning solution may be used to clean the external surfaces and internal connecting surfaces of the attachments and tubes. Straight attachments, footed attachments and telescoping straight tubes have special cleaning brushes sized to the attachment's or telescoping tube's internal diameter. Push the brush wet with surgical instrument cleaning solution through the attachment or telescoping tube from rear to front to loosen and remove debris trapped inside. - Move any moveable parts back and forth to allow solution to thoroughly clean attachment, e.g., sleeve on footed attachment, perforator attachment. - Rinse thoroughly with tap water. - Thoroughly dry attachments. An air gun may be used to blow moisture out from rear to front of attachment. - Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: - Holding the can approximately 10-15 cm (3-6 in.) away from the attachment, spray all components that move, rotate, or slide with three quick squirts. - Articulate movable components to ensure proper lubrication. - Remove excess lubricant with a clean cloth.
Disinfection	No particular requirements
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap

Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization:				
	Cycle:	Gravity	Pre-vac	Pre-vac (FR/WHO)*	Pre-vac (UK)*
	Temperature:	132°C	132°C	134°C	134°C
	Time:	25 min.	4 min.	18 min.	3 min.
	Drying:	15 minutes	15 minutes	20 minutes	10 minutes
*Items contaminated with Transmissible Spongiform Encephalopathies (TSE) agents may be decontaminated using steam autoclave at a temperature of 134–137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles as referenced in NHS Estates HTM 2010 parts 4 & 6: Appendix 2, Items Contaminated With TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies. Medtronic recommends that all products used on a patient suspected or confirmed with a TSE diagnosis be incinerated.					
STERRAD Sterilization: Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.					
100% EtO Sterilization Parameters: Preconditioning: 51–59°C, 70 ±5% relative humidity, 30 min. Temperature: 51–59°C Relative Humidity: 70 ±5% Ethylene oxide concentration: 725 ± 25mg/L Gas exposure time (full-cycle): 4 hours Aeration: 18 hours at 51–59°C					
Steris: Do not use liquid peracetic acid sterilization due to immersion procedure.					
Maintenance, Inspection and Testing	Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. Verify functionality prior to re-use.				
Storage	Store with other sterile devices.				

NOTE

The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

Legend Pneumatic Control Unit / Instrument Case

Warnings and Precautions	Do not soak/submerge Legend devices. Do not use ultrasound to clean Legend devices. Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. Do not sterilize the Legend pneumatic control unit.
Limitations	Verify functionality prior to re-use.

INSTRUCTIONS

Point of Use	No particular requirements.
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.
Preparation for decontamination	No particular requirements.
Cleaning: Automated (DO NOT use ultrasonic washer)	Not validated.
Cleaning: Manual	Wipe the pneumatic control unit and instrument case with a cloth dampened with surgical instrument cleaning solution after each use.
Disinfection	No particular requirements. Follow hospital procedures.

Packaging	No particular requirements. Follow hospital procedures.
Sterilization	Do not sterilize pneumatic control unit.
Maintenance, Inspection and Testing	Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. Verify functionality prior to re-use.
Storage	Do not store with sterile devices.

NOTE

The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

Transmissible Spongiform Encephalopathies (TSE) Return Policy

Medtronic Powered Surgical Solutions (MPSS) will not authorize or accept the return of MPSS products that directly contact patients or is contaminated with a patient's body fluids who is suspected or confirmed with a Transmissible Spongiform Encephalopathies / Creutzfeldt-Jakob Disease (TSE/ CJD) diagnosis. Furthermore, MPSS recommends that all MPSS products used on a patient confirmed with a TSE diagnosis be incinerated. Contact your Sales Representative for replacement of product incinerated under this policy or for temporary equipment while original equipment is quarantined. Contact MPSS Regulatory Affairs Department for additional information regarding TSE contamination. MPSS dissecting tools used on a patient suspected of a TSE/CJD diagnosis must be incinerated.

If TSE/CJD is excluded as a diagnosis, the quarantined reusable equipment may be returned for use after appropriate cleaning, decontamination and sterilization. Hospital personnel should contact their infection control personnel for current procedures and policy for reusable equipment processing when suspect of contamination with Creutzfeldt-Jakob Disease (CJD) or other Transmissible Spongiform Encephalopathies (TSE).

Reusable devices that have been used on patients with suspected Creutzfeldt-Jakob Disease (CJD) or other Transmissible Spongiform Encephalopathies (TSE) should be quarantined and not reused until diagnosis is confirmed or excluded. Reusable devices should be quarantined after having been cleaned, decontaminated, sterilized and packed in a ridged sealed container until final diagnosis.

STRAIGHT ATTACHMENTS

To Install:


1. Slide straight attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when attachment is fully seated.

NOTE

A Legend Attachment will not seat on a Legend Motor if the markers on the collet flats are not in alignment.

2. Insert dissecting tool into attachment with a slight rotational motion. A tactile and audible click is observed indicating that the tool is fully seated.

CAUTION

Match the nomenclature and color code on the Legend dissecting tool packaging to the same nomenclature and color code on the Legend attachment.

3. Rotate attachment to the "☑" position on the motor case. Gently pull on dissecting tool shaft to ensure proper installation.

NOTE

The Legend Motors will not run properly unless the attachment is in the "☑" position.


Smoke may be generated if attachment is not in the locked "☑" position.

To Remove:

1. Hold motor in palm of hand. Rotate attachment to the "☐" position.

2. Remove dissecting tool from the attachment and discard dissecting tool.

3. Use thumb and index finger to lift attachment off motor.

ANGLED ATTACHMENTS

To Install:

1. Slide angled attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when attachment is fully seated.
2. Rotate attachment to "🔒" position on motor case.
3. With "Tool Lock" in the unlocked position, insert dissecting tool into angled attachment with a slight rotational motion. A tactile and audible click is observed indicating that the tool is fully seated.
4. Rotate Tool Lock to the "🔒" position. Gently pull on dissecting tool shaft to ensure proper installation.

NOTE

A dissecting tool may be installed and locked in the attachment before the angled attachment is installed onto the motor.



Angled and straight attachments with the same length, marking and color band share the same dissecting tools.

The Legend Motors will not run properly unless the attachment is in the "🔒" position.

Smoke may be generated if attachment is not in the locked "🔒" position.

To Remove:

Rotate the Tool Lock to the "🔓" position to remove the dissecting tool from the attachment. Rotate the attachment to the "🔓" position and lift attachment off of the motor.

⚠ CAUTION

When using an angled attachment, hold the handpiece assembly by the attachment so that the attachment does not inadvertently loosen from the handpiece.



ANGLED DOUBLE LOCK ATTACHMENTS

To Install:

1. Slide the attachment over the motor collet, aligning the triangular markers.
2. Pull the attachment towards the motor and turn the attachment to the "☰" position.
3. With the "tool lock" in the "☐" position, insert the dissecting tool into the attachment with a slight rotational motion. You will feel the tool click into place when it is fully seated.
4. Turn the tool lock to the "☒" position. Gently pull on the shaft of the dissecting tool to verify proper installation.
5. On Variable Exposure Attachments, use the TUBE adjustment ring to adjust the exposure of the dissecting tool. With the tool pointing away from you, turn the ring to the right to increase the length of the tube, thereby decreasing the exposure of the tool. Turn the ring to the left to decrease the length of the tube, thereby increasing the exposure of the tool.

To Remove:

1. To remove the attachment, hold the motor in the palm of your hand, and push the sleeve on the attachment distally while turning the attachment to the "☒" position.
2. Release the sleeve and remove the attachment.

CURVED BUR ATTACHMENTS

⚠ WARNING

Use adequate irrigation and keep the cooling sleeve soaked during dissection. Inadequate irrigation may cause thermal necrosis.

Do not modify the bur. Bending or prying may break the bur, causing harm to the patient or operating room staff.

Excessive pressure applied to the bur may cause bur damage. If this occurs, use extreme care to ensure that all fragments of the bur are removed from the patient.

Disposable devices are for single-use only. Do not attempt to sterilize disposable devices. The disposables are packed sterile and are not intended for repeat use. To prevent contamination, use only once.

Test for bur wobble (eccentricity) at the desired speed prior to use. Select a new bur or reduce speed if wobble is observed prior to use or during the procedure. Bur wobble may cause patient injury.

⚠ CAUTION

Curved burs are not designed for variable tool exposure. Do not attempt to adjust tool exposure, as this may damage the device.

When operating or testing the motor, ensure the bur is properly inserted and locked into AT10 attachment.

Installation

1. Turn the tube locking ring and the tool locking ring on the AT10 (V01) Telescoping Attachment into the “” position.
2. Slide the curved bur into the AT10 (V01), until the hub is fully seated.
3. Gently press on the head of the dissecting tool to lock it into place. You will hear a click when the curved bur is fully seated.
4. Turn both locking rings to the “” position on the AT10 (V01).
5. Gently pull on the curved bur to verify proper installation.
6. Refer to your motor’s *Instruction Manual* for information on setting up the AT10 (V01) attachment and motor.
7. Briefly run the motor with the curved bur installed, checking for bur wobble or excessive vibration. If either occurs, perform the corrective actions described in the *Troubleshooting* section of these instructions.

Cooling the Bur

1. Prior to initial use, soak the cooling sleeve by dipping it into a cup of saline or DI water, as shown below.

Wetting of the cooling sleeve prior to cutting

2. During use, maintain copious irrigation of the cooling sleeve and bur tip by dribbling saline or DI water along the entire length of the cooling sleeve and bur tip.

Operation

1. Activate the motor and gently press the curved bur against the bone to begin dissection.
2. Use a light sweeping motion to continue removing bone.

Removal

1. Turn both AT10 (V01) locking rings to the “” position and pull the bur out of the attachment.
2. Discard the curved bur, according to hospital procedures.

Reuse and Cleaning

Do not reuse.

Sterility

Each bur is gamma-sterilized and is not intended for repeated use. Do not attempt to re-sterilize.

Cutting Time

Continuous cutting for extended periods may cause the device to heat to an uncomfortable temperature. To avoid this, limit cutting to the recommended times below:

	Maximum Continuous Cutting Time	Maximum Total Cutting Time
1 & 2 mm burs	2 minutes	12 minutes
3 & 4 mm burs	3 minutes	12 minutes

VARIABLE EXPOSURE ATTACHMENTS

⚠ CAUTION

The Legend Variable Exposure attachments can be distinguished from standard attachments by the dual color bands on the attachment. Match the color band on the attachment to the color code on the dissecting tool packaging.

⚠ WARNING

Surgeons should familiarize themselves with the performance of dissecting tools before use, and should explore the effect of various levels of tool exposure on dissection stability. If the tool exhibits excessive chatter, vibration, or movement, decrease the tool exposure.

Dissecting tool size and geometry may contribute to excessive vibration at certain speeds. Increase or decrease speed by adjusting the foot/finger control, or by changing the operating pressure or console speed setting. If necessary, use a different dissecting tool.

To Install:

1. Refer to the appropriate section of this manual on installing fixed or angled attachments.
2. After installation, use the TUBE adjustment ring to adjust the exposure of the dissecting tool. With the tool pointing away from you, turn the ring to the right to increase the length of the tube, thereby decreasing the exposure of the tool. Turn the ring to the left to decrease the length of the tube, thereby increasing the exposure of the tool.

⚠ WARNING

Do not use the Variable Exposure Attachment if the TUBE adjustment ring spins freely or fails to click into place with each adjustment, as the exposure may change without warning.

⚠ CAUTION


Make sure that the tool lock (angled attachments only) and the attachment lock are still in the "🔒" position after each adjustment of the tool exposure. Attempting to increase the exposure too far may result in the attachment becoming unlocked. Accidentally turning the tool lock may result in reduced speed and/or overheating of the attachment.

⚠ WARNING

Do not use the end of the tube as a depth gauge or depth stop.

To Remove:

1. Removal is the reverse of installation.
2. When cleaning, clean the attachment completely, first without adjusting the tube length, then with the tube fully extended, and with the tube fully retracted.

ROTATING FOOTED ATTACHMENTS

NOTE

Rotating and fixed footed attachments with the same length, marking, and color band share the same dissecting tools.

To Install:

1. Insert the dissecting tool into the motor collet with a slight rotational motion. You will feel and hear the tool click into place when it is fully seated.
2. Slide the attachment over the dissecting tool and onto the motor, aligning the triangular markers. You will feel and hear the attachment click into place when it is fully seated.
3. Turn the attachment to the “” position.

NOTE

The footed end of the attachment now has 360° of unrestricted rotation.

To Remove:

Removal is the reverse of installation.

FOOTED ATTACHMENTS



Footed Attachments:	Use only the following Dissecting Tools:	Reference to Classic:
F1-B5	F1/8TA15 F1/8TA15S	B- 5 22S2CR (Mednext)
F2-B1	F2/8TA23 F2/8TA23S	B-1 B-1SP
F3-S	F3/9TA30	S-1

⚠ WARNING

DO NOT use a Legend Attachment if any part of the attachment appears to be bent, loose, missing or damaged. Excessive pressure or improper handling, such as bending or prying, with the attachment or dissecting tool may cause injury to patient, operator and/or operating room staff. If a Legend Attachment requires servicing or refurbishing, return attachment to MPSS.

To Install:

1. Insert dissecting tool into motor collet with a slight rotational motion. A tactile and audible click is observed indicating that the dissecting tool is fully seated.
2. Slide footed attachment over dissecting tool onto motor aligning triangular markers on attachment and motor case.

3. Pull footed attachment towards motor and rotate attachment to the "■" position on motor case.


NOTE

Legend footed attachments have a unique feature of actuating a distal moving sleeve to remove the attachment. This will prevent accidental disengagement of attachment from motor during use.

The Legend Motors will not run properly unless the attachment is in the "■" position.

To Remove:

1. To remove Legend footed attachment, hold motor in palm of hand. Push sleeve on footed attachment distally while rotating attachment to the "■" position on the motor case. Release sleeve.
2. To avoid injury from the dissecting tool, use thumb and index finger to cautiously and slowly lift attachment off of motor and away from the dissecting tool.
3. Pull dissecting tool out of motor collet and discard tool.

⚠ WARNING

Remove Legend Footed Attachments cautiously and slowly to avoid injury to the operator.

CONTRA-ANGLE ATTACHMENT 16-MF


Dynamic Operating Pressure: 80 psi (5.5 bar)

CAUTION

DO NOT run the 16-MF attachment with operating pressure above 80 psi (5.5 bar). This may cause over-heating and damage to internal gears of attachment.

To Install:


- 1. Decrease pressure at the compressed gas source to 80 psi (5.5 bar).
Adjust pressure as needed by lowering pressure at gas source, then push down pressure relief on Pneumatic Control Unit to exhaust excess pressure in hoses. Verify at supply pressure.
- 2. Slide Contra-Angle Attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when attachment is fully seated.
- 3. Rotate attachment to the “” position on the motor case.
- 4. Rotate attachment head's lever laterally to open position. Insert dissecting tool and return lever to closed position. Gently pull on dissecting tool shaft to ensure proper installation.

Lateral View

Lever Open

Lever Closed

NOTE

A dissecting tool may be installed and locked in the attachment before the Contra-Angle attachment is installed on the motor.

The Legend Motors will not run properly unless the attachment is in the “” position.

To Remove:

Rotate lever on attachment head laterally to remove the dissecting tool. Discard dissecting tool. Rotate the 16-MF attachment to the “” position and lift attachment off of the motor.

METAL CUTTING ATTACHMENT

The Metal Cutting Attachment uses the tungsten carbide or diamond wheel dissecting tools. All metal cutting dissecting tools have an "MC" attachment prefix in their nomenclature, e.g., MC254, MC30.

Metal cutting dissecting tools cannot be installed into any other Legend Attachment.

WARNING

DO NOT use metal cutting dissecting tools on bone.

For metal transection, observe the following Safety Precautions:

- Eyewear protection is essential.
- Irrigate copiously to cool the cutting surfaces.
- Protect the wound site from metal debris.
- Use clamp or grasping device to control loose fragments during transection of any metal component.

To Install:

1. Slide Metal Cutting Attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when attachment is fully seated. Rotate attachment to the "🔒" position on the motor case.

2. With tool lock unscrewed several turns and holding the attachment UPRIGHT, insert dissecting tool into attachment. Rotate dissecting tool until it drops into position and is fully seated. A tactile click is observed indicating that the tool is fully seated.

3. Firmly and completely tighten the tool lock. Pull on the tool shaft to ensure that it is properly installed.

NOTE

A dissecting tool may be installed and locked in the metal cutting attachment before the attachment is installed onto the motor.

The Legend Motor will not run properly unless the attachment is in the "🔒" position.

To Remove:

1. Unscrew the tool lock with several turns, then withdraw the dissecting tool.
2. Rotate attachment to the "🔓" position on the motor case and lift attachment off of the motor.

TELESCOPING ATTACHMENT

To Install:


1. Slide Telescoping Attachment over motor collet aligning triangular arrows on attachment and motor case. A tactile and audible click will be observed when attachment is fully seated.

Rotate attachment to the “” position on motor case.

2. Rotate TUBE Locking Ring towards “” by turning counterclockwise. Insert base end of selected telescoping tube into attachment. To lock tube in place, turn the TUBE Locking Ring clockwise until finger tight.

TUBE
Locking
Ring



NOTE

DO NOT over tighten.

3. Be sure that the TOOL Locking Ring is in the “” position. Insert the dissecting tool into the telescoping tube. Tactile feedback is felt when the dissecting tool is fully seated. Turn the TOOL Locking Ring to the “” position. Verify that the tool is in place by gently pulling on the tool.
4. If tube position needs to be changed, rotate TUBE Locking Ring towards “”, re-position tube, then rotate TUBE Locking Ring towards “”. Gently pull on dissecting tool, then the tube, to ensure proper installation.

TOOL
Locking
Ring

NOTE

A dissecting tool may be installed and locked in the attachment before the telescoping attachment is installed on the motor.


Curved and straight telescoping tubes with the same length, marking and color band share the same dissecting tools.

Telescoping tubes are disposable following multiple uses and should be discarded when heat or excessive vibration is noticed, or when insertion of the tools becomes difficult.

The Legend Motor will not run properly unless the attachment is in the “” position.

To Remove:

1. Rotate the TUBE Locking Ring towards “”.
2. Rotate the TOOL Locking Ring to the “” position.
3. Pull telescoping tube out of attachment.
4. Rotate the attachment to the “” position on the motor case and lift attachment off of motor.

WARNING

Dissecting tool flutes are sharp and may perforate surgical gloves. Dissecting tools for the Telescoping Attachment may be grasped with a hemostat to aid in installation and removal.

PERFORATOR ATTACHMENT

The Perforator Attachment has a Hudson chuck to drive any device with a Hudson shank, i.e., cranial perforator device.

To Install:

1. Slide Perforator Attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when Perforator Attachment is fully seated.

2. Rotate Perforator Attachment to the "L" position on the motor case.

3. To install a cranial perforator device with a Hudson shank, pull back proximally on the collar of the Perforator Attachment. Insert device and release collar to its original position.

NOTE

Cranial perforator device may be installed in the attachment before the perforator attachment is installed on the motor.

MPSS does not provide cranial perforator devices.

The Legend Motors will not run properly unless the attachment is in the "L" position.

To Remove:

Pull back proximally on the collar of the Perforator Attachment to remove cranial perforator device. Rotate the Perforator Attachment to the "L" position and lift attachment off of the motor.

WARNING

CONSULT the cranial perforator device labeling for the recommended speed specifications.

Maximum Speed of Perforator Attachments with Legend Pneumatic and EHS Motors

Pneumatic Motors		
Gas Pressure (Dynamic)	AD01 Output Speed (Max.)	AD03 Outut Speed (Max.)
80 psi	660 rpm	850 rpm
100 psi	820 rpm	1050 rpm
120 psi	890 rpm	1140 rpm
EHS Motor		
EHS Console Setting	AD01 Output Speed (Max.)	AD03 Output Speed (Max.)
60,000 rpm	645 rpm	830 rpm
70,000 rpm	745 rpm	965 rpm
72,000 rpm	770 rpm	995 rpm
74,000 rpm	790 rpm	1020 rpm
75,000 rpm	805 rpm	1035 rpm

JACOBS® CHUCK ATTACHMENT

The Jacobs Chuck Attachment is a non-cannulated 5/32" chuck with key for drilling.

To Install:

1. Slide Jacobs Chuck Attachment over motor collet aligning triangular markers on attachment and motor case. A tactile and audible click will be observed when Jacobs Chuck Attachment is fully seated.

2. Rotate Jacobs Chuck Attachment to the "🔒" position on the motor case.

3. To install a drill bit, turn Jacobs key to open ridged collar. Insert drill bit and tighten collar with key.

NOTE

A drill bit may be installed in the attachment before the Jacobs Chuck attachment is installed on the motor.

MPSS does not provide drill bits.

The Legend Motors will not run properly unless the attachment is in the "🔒" position.

To Remove:

Use Jacobs key to open collar. Remove and discard drill bit in an appropriate container. Rotate the Jacobs Chuck Attachment to the "🔓" position and lift attachment off of the motor.

BONE MILL ATTACHMENT

The Legend Bone Mill Attachment is composed of a non-sterile base and gamma-sterilized, single use disposable components. The gamma-sterilized single use disposable components consist of one bowl, one cap (with LOW and HIGH settings), and one spatula. A Legend or MR7 motor provides the power to drive the Bone Mill Attachment. The Bone Mill Attachment is intended to mill bone into particles 1 to 5mm in size.

WARNING


DO NOT operate the Bone Mill Attachment without the bowl, and cap secured in place.

DO NOT come in contact with the interior of the disposable bowl and cap during bone milling.

DO NOT use the Bone Mill disposable components for more than one surgical procedure, as this may cause cross-contamination and affect patient safety.

Use only Medtronic Legend/MR7 and Legend Bone Mill disposables in combination with the Bone Mill Attachment. Use of other devices will void the manufacturer's warranty.

Legend Bone Mill Attachment Components

Using the Bone Mill Attachment

1. Place the bowl onto Bone Mill base and secure the latches.
2. Lay harvested bone specimens into the bowl.

CAUTION

The addition of any liquid to the milling process should be administered in small volumes to prevent leakage.

DO NOT use Bone Mill disposables to process more than 6 bone pieces approximately 1.5 cc in size at one time.

DO NOT use Bone Mill disposables to process more than a total of 40 cc of bone during any one surgical procedure.

3. Secure the cap onto the bowl by rotating counter-clockwise. Do not activate the Bone Mill Attachment without securing the cap.
4. Use the HIGH cap position when starting with large pieces of bone or large batch amounts. Use the LOW cap position to finish milling a batch or when starting with small pieces or batch amounts.

5. Insert the motor into the motor connection on the Bone Mill base and rotate the motor to the "■" position.

CAUTION

DO NOT run the Legend Bone Mill Attachment with a pneumatic motor at an operating pressure below 80 psi (5.5 bar) or above 120 psi (8.3 bar).

6. Activate the motor to mill bone. Continue to run the Bone Mill Attachment until desired bone particle size is achieved. Adjust motor speed to obtain desired particle size.

NOTE

If bone does not process, deactivate the motor, rotate the motor to the "■" position, and remove the cap. Reinsert the cap in the HIGH position and activate the motor. If bone spins without milling or jams, divide bone chunks into two batches and process separately.

7. Process the bone for 5-15 seconds and visually check to gauge quality of milling before continuing.
8. Continue to process bone in 5-15 second batches until the majority of particles reach desired size.
9. To remove bone particles stop operation of the motor and rotate it to the "■" position. Remove the cap by rotating it clockwise. Use the spatula inside the bowl and cap to remove bone particles.

NOTE

Milling time is dependent on the density of the patient's bone, the ratio of cancellous and cortical bone, the size of the individual particles and the batch amount of the bone to be processed. Always check bone after 10-15 seconds to gauge quality of milling.

NOTE

If the blade becomes inoperable, stop operation of the motor and rotate the motor to the “” position. Remove the cap and use the spatula to free bone and/or blades.

After removing bone particles from the bowl, the milling process may be repeated with additional bone specimens, as needed. (Up to a total of 40 cc of bone may be processed by the Bone Mill Disposable during any one surgical procedure.)

After Use

- 1. Rotate the motor to the “” position.
- 2. Unlatch the Bone Mill disposables from the base.
- 3. Discard the bowl, cap, and spatula per hospital procedure.
- 4. Clean the Bone Mill base according to the directions below.
- 5. Clean the motor according to its *Instruction Manual*.
- 2. Bone Mill bowl, cap, and spatula are single use devices. Opened but unused Bone Mill disposables may not be re-sterilized under any sterilization method.

Preventative Maintenance

The Bone Mill Attachment should be returned for factory maintenance annually (approximately 100 procedures).

If the Bone Mill Attachment does not function according to these instructions, call your Medtronic Regional Distributor or Medtronic Neurologic Sales Representative. If necessary, return the Bone Mill Attachment base to your Distributor or Sales Representative to be repaired.

Contact your Medtronic Neurologic Technologies sales representative, or call:

Medtronic Powered Surgical Solutions Technical Services : 8 a.m. to 6 p.m. (Central Standard Time)
Telephone: 800-643-2773 or 817-788-6400

Outside the U.S.A., contact your nearest Medtronic Regional Distributor or Medtronic Neurological Sales Representative or call U.S.A. Telephone 817-788-6400

NO LATEX POLICY:

MPSS verifies that all Legend equipment and products identified in this manual manufactured by and/or for MPSS do not contain latex. MPSS packaging materials, labels, pouches, package inserts or similar items do not contain latex.

Material Safety Data Sheet:

Material Safety Data Sheet on lubricant used in Legend Lubricant/Diffuser is available by contacting Regulatory Affairs Department or Technical Service at telephone numbers listed above.

LEGEND SYSTEM TROUBLESHOOTING

If Legend components require servicing or refurbishing, return components to MPSS for quality assured service by factory-trained personnel who will utilize genuine MEDTRONIC MIDAS REX LEGEND parts as required.

All Legend System components returned for servicing or refurbishing should be properly cleaned and sterilized prior to shipping.

PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
LEGEND MOTORS		
Motor does not run or is low on power.	Hoses not properly connected.	Check all connections from compressed gas source to motor.
	Operating pressure inadequate.	Check gas supply to Pneumatic Control Unit at Supply Pressure gauge. Pressure may need to be increased according to compressed gas requirements.
	Attachment not properly installed and locked onto the motor (possible smoking condition).	Remove and re-install attachment and dissecting tool to ensure proper installation and locking of attachment onto motor.
	Motor sterilized with single use, disposable Lubricant/Diffuser cartridge attached to the motors lubricant/diffuser housing.	DO NOT use. Return motor to MPSS to be refurbished.
	Foot Pedal on Pneumatic Control Unit not properly functioning.	Check for obstruction under foot pedal. If foot pedal continues not to function properly, return the Pneumatic Control Unit to MPSS to be refurbished.
	Motor stall or will not spin.	Manually spin the dissecting tool then activate the motor. If motor continues to stall, then return motor to MPSS to be refurbished.

LEGEND SYSTEM TROUBLESHOOTING (continued)

PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
LEGEND MOTORS (continued)		
 Legend Gold Touch will not run after checking possible causes listed previously for all Legend motors.	Finger/Foot control on the Pneumatic Control Unit not activated.	Push finger/foot control on Pneumatic Control Unit to the “  ” position.
	Safety slide in “○” position.	Move safety slide to “ ” position.
Motor continues to run.	Pneumatic Control Unit is not properly functioning.	Depressurize system. Return Pneumatic Control Unit to MPSS to be refurbished.
	Finger control is not properly functioning	Return motor to MPSS to be refurbished.
Motor is missing a component or a component has detached from the motor.	Motor failed.	DO NOT use. Return motor to MPSS to be refurbished.
Motor rotates, but an abnormal noise is heard.	Lubrication inadequate.	Check for proper installation of the Legend Lubricant/Diffuser. If problem persists, replace Legend Lubricant/Diffuser. If problem continues, return motor to MPSS to be refurbished.
	Motor’s exhaust hose damaged or cut, or internal pressure hose detached.	DO NOT use. Depressurize system. Return motor to MPSS to be refurbished.
	Bearings of motor or attachment worn.	Change attachment to isolate location of problem to the attachment or motor. Return motor and/or attachment to MPSS to be refurbished.
	Attachment not in proper locked position (possible smoking condition).	See instructions for proper locking of attachments.
Motor/collet displays excessive rust.	Incorrect cleaning or sterilization method.	Return to MPSS for refurbishing/purchase.
	Use of chlorine based or corrosive cleaning agents.	
	Aging.	

LEGEND SYSTEM TROUBLESHOOTING *(continued)*

PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
LEGEND MOTORS <i>(continued)</i>		
Motor is heating to an uncomfortable temperature to touch/hold.	Inadequate "cool down" period following sterilization.	Motor must be allowed to cool down following steam sterilization. Run motor for 2 minutes. If motor continues to overheat, return the motor to MPSS to be refurbished.
	Lubrication inadequate.	Check for proper installation of a new Legend Lubricant/Diffuser into motor's lubricant/diffuser housing.
	Attachment transferring heat to motor.	Change attachment to isolate location of heating to the attachment or the motor. If attachment is source of heat, return attachment to MPSS to be refurbished. If motor is source of heat, return the motor to MPSS to be refurbished.
	Heavy side loading during dissection.	Discontinue use and rest the motor by using intermittently or wrap the motor interface with a moist sterile towel. If motor continues to overheat, return the motor to MPSS to be refurbished.
	Inadequate irrigation.	Ensure adequate irrigation to surgical site during bone dissection.
Motor is leaking lubricant.	Motor failed.	DO NOT use. Return motor to MPSS to be refurbished. If lubricant entered surgical site, treat patient per physician's discretion.
Oily film on green exhaust hose of motor following sterilization.	Oily film on hose due to pressure and/or temperature differentials following sterilization.	Wipe the external surface of the motor hose with a sterile cloth prior to use. After wiping with sterile cloth, if motor continues to have oil on the hose, DO NOT use. Return the motor to MPSS to be refurbished.
Hose explodes.	Hose pinched, kinked, obstructed or stepped on.	DO NOT use. Return motor to MPSS to be refurbished/purchase.

LEGEND SYSTEM TROUBLESHOOTING (continued)

PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
LEGEND ATTACHMENTS OR TELESCOPING TUBES		
Attachment or Telescoping Tube has uncomfortable temperature to touch/hold.	Heat from worn attachment/tube bearings.	DO NOT use. Try another attachment/tube. Return over-heated attachment to MPSS to be refurbished. Telescoping Tubes are multi-use disposable. If problem is resolved with a new Telescoping Tube, discard the over-heated tube.
	Attachment/tube unclean due to improper cleaning procedures.	Check that appropriate cleaning procedures are being followed.
	Heavy side loading during dissection.	Discontinue use and rest the attachment by using intermittently, try another identical attachment or wrap the attachment interface with a moist sterile towel. If attachment continues to overheat, return the attachment to MPSS to be refurbished/purchase.
Color band on Attachments/Telescoping Tubes fades or discolors.	Incorrect cleaning or sterilization method.	Use nomenclature markings on the attachment to match with a corresponding dissecting tool or return to MPSS for refurbishing/purchase. Telescoping Tubes are multi-use disposable.
	Use of chlorine based or corrosive cleaning agents.	
	Aging.	
Attachment/Telescoping Tube displays excessive rust.	Incorrect cleaning or sterilization method.	DO NOT use. Return attachment to MPSS to be refurbished/purchase. Dispose of telescoping tube. Telescoping Tubes are multi-use disposable.
	Use of chlorine based or corrosive cleaning agents.	
	Aging.	
Attachment/Telescoping Tube is bent, loose, damaged or missing a component.	Attachment mishandled, failed due to extended use or excessive force applied during use.	DO NOT use. Return attachment to MPSS to be refurbished/purchase. Dispose of telescoping tube. Telescoping Tubes are multi-use disposable.

LEGEND SYSTEM TROUBLESHOOTING *(continued)*

PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
LEGEND ATTACHMENTS OR TELESCOPING TUBES <i>(continued)</i>		
Attachment will not properly seat on motor.	Motor collet flats are not in proper alignment.	Use Legend Motor Wrench to rotate the flat closest to the motor case until its marker is aligned with the marker on the flat farthest away from the motor case.
Attachment has excess lubrication.	Over lubrication during cleaning process.	Visually inspect and wipe excess lubrication.
Difficulty removing tool from attachment.	Aging of attachment. Improper cleaning. Use of re-processed tools. Use of unauthorized re-furbisher.	Return to MPSS for refurbishing/purchase.
Footed attachment has a component missing from leg/foot area.	Attachment damaged by dissecting tool drilling out part or all of leg/foot area.	DO NOT use. Return attachment and motor to MPSS to be refurbished/purchase.
Overheating of 16-MF Contra-Angle Attachment.	Contra-Angle attachment operates by a set of internal gears to engage the drive shaft. It is normal for some heat to be generated approximately 2 cm from the distal end of the attachment and at the right angle head.	Verify pressure setting is at 80 psi (5.5 bar). DO NOT use if heat continues or is excessive. Return attachment to MPSS to be refurbished/purchase
Perforator running too slow.	Pressure set incorrectly.	Check pressure setting at foot control.
LEGEND DISSECTING TOOLS		
Dissecting Tool flails.	Device other than a Legend Dissecting Tool used.	DO NOT use. Replace with a Legend Dissecting Tool.
	Worn attachment/tube bearings.	DO NOT use. Try another attachment/tube to isolate location of problem. If flail is produced by a specific attachment, return attachment to MPSS to be refurbished/purchase. Dispose of telescoping tube. Telescoping Tubes are multi-use disposable.
	Attachment/tube and tool combination improper.	DO NOT use. Check for the correct nomenclature and color code on the dissecting tool packaging to the same nomenclature and color code on the attachment/tube.

LEGEND SYSTEM TROUBLESHOOTING (continued)

PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
LEGEND DISSECTING TOOLS (continued)		
Dissecting Tool flails. (continued)	If proper attachment/tube and tool combination is being used, motor may be damaged. Dissecting tool's size and geometry may contribute to flailing at certain speeds.	DO NOT use. Return motor to MPSS to be refurbished. Increase or decrease speed by changing pressure setting or foot/ finger control. DO NOT use if flailing persists. Change to a different dissecting tool.
Excessive vibration.	Dissecting tool's size and geometry may create excessive vibration at certain speeds.	Increase or decrease speed by changing pressure setting or foot/finger control. Change to a different dissecting tool.
Dissecting tool sterile packaging damaged. Tool outside the tube or tube opened.	Shipping damage - Dissecting tool is no longer sterile.	DO NOT use, contact MPSS Customer Service or Sales representative.
Dissecting tool sterile packaging opened, tool not used or contaminated.	Dissecting Tool not used during surgical procedure - Dissecting tool is no longer sterile.	See instruction for re-processing tools under sterilizations guideline section.
Expired Shelf Life.	Old product - Dissecting tool may no longer be sterile.	DO NOT use, contact MPSS Customer Service or Sales representative.
Dissecting tool will not seat properly in motor collet or angled attachment collet.	Debris in collet of angled attachment or motor. Device other than Legend Dissecting Tool being used.	Clean angled attachment or motor thoroughly as per instructions. If cleaning does not correct the problem, DO NOT use. Return angled attachment or motor to MPSS to be refurbished. DO NOT use. Replace with a Legend Dissecting Tool.
Dissecting tool rusts.	Incorrect storage. Reprocessed unused dissecting tool.	DO NOT use, return to MPSS. DO NOT use.
Dissecting tool fractures or breaks.	Excessive pressure. Reprocessed tool. Tool failure.	Remove all tool debris from patient. Check if reprocessed tool was used. Contact MPSS Customer Service or Sales Representative.
Dissecting tool dull.	Extended use. Reprocessed tool. Incorrect geometry.	Change to a new tool. Check if reprocessed tool was used. Contact MPSS Customer Service or Sales Representative.

REFURBISHING OR REPAIRS

When the Midas Rex Legend Pneumatic High Speed System requires servicing or refurbishing, contact Medtronic Powered Surgical Solutions Repair Services for a return authorization and instructions for returning the equipment to Medtronic Powered Surgical Solutions, Fort Worth, Texas. Medtronic Powered Surgical Solutions provides quality assured service by factory-trained personnel who will utilize genuine Midas Rex Legend parts as required.

Contact:

Medtronic Powered Surgical Solutions Repair Services:

(800) 335-9557 or (817) 788-6440

mmrcustomerservice@medtronic.com

Peak performance, reliability and maximum service life from your Midas Rex Legend System may be assured by using only those Midas Rex Legend products for your Midas Rex Legend System that are manufactured by and sold through Medtronic Powered Surgical Solutions, Fort Worth, Texas. While Medtronic Powered Surgical Solutions guarantees complete compatibility among its products within a specific product line, the dissecting tools are designed for single-use only, and Medtronic disclaims any responsibility when reprocessed dissecting tools are used. If you would like more information about the patient and product risks associated with reprocessed tools, please contact the number or e-mail address listed above.

Due to safety and environmental concerns, Medtronic Powered Surgical Solutions requests the return of pneumatic high speed motors for proper disposal at the end of the product life cycles.

LEGEND PNEUMATIC HIGH SPEED PREVENTIVE MAINTENANCE/SERVICE PROGRAM

The Midas Rex Legend Pneumatic High Speed System Preventive Maintenance Manual has been developed to assist you in getting the greatest ownership value from your Midas Rex Legend Pneumatic High Speed System, while helping to maximize its performance, safety and reliability. The scheduled preventive maintenance/service program is in addition to the required routine cleaning after each use. Please refer to the preventive maintenance manual for the specific steps necessary to maintain Midas Rex Legend Pneumatic High Speed System.

MANUFACTURER'S RECOMMENDATION

It is strongly recommended that surgeons, operating room personnel, and service personnel become thoroughly familiar with all Midas Rex Legend Systems prior to use in surgery. To assist in this effort, Medtronic Powered Surgical Solutions offers several Hands-On Workshops. Contact your Medtronic Neurologic Technologies Sales Representative or your Medtronic Regional Distributor for schedules.

MIDAS REX LEGEND PNEUMATIC HIGH SPEED SYSTEM LIMITED WARRANTY*

(U.S. Customers Only)

- A. This Limited Warranty provides the following assurance to the purchaser of a Medtronic Midas Rex Legend Pneumatic High Speed System. The Midas Rex Legend Pneumatic High Speed System includes the motor, foot control, instrumentation cases and trays (hereafter referred to as Motor Components), the attachments, dissecting tools and other accessories not listed above and telescoping tubes (hereafter referred to as Accessory Components) and jointly referred to as Legend Pneumatic High Speed System, unless specifically noted.
- (1) Should a Motor Component fail to function within normal tolerances due to a defect in materials or workmanship within a period of one year from the date of shipment of the new Motor Component, or 180 days from the date of shipment of the serviced or refurbished Motor Component to the purchaser (but not less than the one year period from the date of shipment of the new Motor Component), or
 - (2) Should an Accessory Component fail to function within normal tolerances due to a defect in materials or workmanship within a period of 90 days from the date of shipment of a new Accessory Component or serviced or refurbished Accessory Component to the purchaser, Medtronic will at its option: (a) issue a credit to the purchaser of the replacement Motor Component or Accessory Component equal to the Purchase Price, as defined in Subsection A(3), against the purchase of any same Motor Component or Accessory Component requested as its replacement, or, (b) provide a functionally comparable replacement Motor Component or Accessory Component at no charge.
 - (3) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original or current functionally comparable, or replacement Motor Component or Accessory Component.
- B. To qualify for this Limited Warranty, these conditions must be met:
- (1) The Legend Pneumatic High Speed System must be used in accordance with the labeling and *Instruction Manual* provided with the products.
 - (2) Only Legend Pneumatic High Speed System dissecting tools that are manufactured, authorized, and quality controlled by Medtronic may be used.
 - (3) No reprocessed dissecting tools may be used.
 - (4) The Legend Pneumatic High Speed System may not be damaged by accident, abuse, misuse or tampering.
 - (5) The Legend Pneumatic High Speed System may not be modified without the written permission of Medtronic.
- C. This Limited Warranty is limited to its express terms. In particular:
- (1) Except as expressly provided by this Limited Warranty, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE LEGEND PNEUMATIC HIGH SPEED SYSTEM TO FUNCTION WITHIN NORMAL TOLERANCES WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, NEGLIGENCE, STRICT LIABILITY, OR OTHER TORT OR OTHERWISE.
 - (2) This Limited Warranty is made only to the purchaser. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM OR OTHERWISE.
 - (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Limited Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid. This Limited Warranty gives the purchaser specific legal rights. The purchaser may also have other rights that vary from state to state.
 - (4) No such express or implied warranty to the purchaser shall extend beyond the period specified in A(1) above.
 - (5) This Limited Warranty shall be the exclusive remedy available to any person.
 - (6) No person has any authority to bind Medtronic to any representation, condition, or warranty, except this Limited Warranty.

*This Limited Warranty is provided by Medtronic Powered Surgical Solutions, 4620 North Beach Street, Fort Worth, Texas 76137-4116. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for terms of the warranty.

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Contact:

Servizi di Riparazione Medtronic Powered Surgical Solutions (MPSS):

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©2014 Medtronic, Inc.
All rights reserved.
Printed in the USA
01/14
175002ML36.RC

Medtronic

Manufactured by:
Medtronic Powered Surgical Solutions
4620 North Beach Street
Fort Worth, TX 76137 USA
Telephone: 817-788-6400
Toll-Free: 1-800-643-2773
Fax: 817-788-6401
www.medtronic.com

Orders:
Telephone: 817-788-6400
Toll-Free: 1-800-433-7080
Fax: 817-788-6489

Repairs:
Telephone: 817-788-6440
Toll-Free: 1-800-335-9557
Fax: 817-788-6441

EC Authorized Representative:
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands

Medtronic Australian Sponsor:
Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde NSW 2113 Australia

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

	Attention, see Instructions for Use.	TUBE	Tube Control
R_x Only	United States federal law restricts this device to sale by or on the order of a physician.		Air Pressure Relief
REF	Reference Number		Use with
LOT	Lot Number	MULTI-USE DISPOSABLE ATTACHMENT	Multi-Use Disposable Attachment
SN	Serial Number	Instrument Case	Instrument Case
	Quantity	Lubricant/Diffuser	Lubricant/Diffuser
	For single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.	Dissecting Tool	Dissecting Tool
≈	Approximately equal to	Attachment	Attachment
STERILE	Non-Sterile	Control Unit	Control Unit
STERILE R	Sterilized by Gamma Irradiation	Refurbished	Refurbished
	Use by date	Accessory	Accessory
	Date of manufacture	REGULATOR	Regulator
	Temperature Limitations	Bone Mill	Bone Mill
	Unlock	MOTOR	Motor
	Lock	Brush	Brush
	On	Adapter	Adapter
○	Off	! USA	USA Only
	Finger-Operated Control		Manufacturer
	Foot-Operated Control		Do not dispose to unsorted municipal waste.
	Compliant with European Council Directive MDD 93/42/EEC.	TOOL	Tool Control
ECT REP	Authorized Representative in the European Community		Consult Instructions for Use

Contents

General Information.....	1
Indications for Use.....	1
Contraindications	1
Special Notices	1
General Safety Precautions.....	2
No Latex Policy	2
System Components.....	3
Non-Disposable Components	3
Disposable Components.....	7
Legend® Dissecting Tools.....	8
General Guidelines For Attachment and Tool Applications.....	9
Setting up the Operating Room	12
Power Source Requirements	12
Pneumatic Connections.....	12
Setting up the MR7 System.....	13
Installing the Oiler Cartridge.....	13
Connecting the Motor	14
Installing an Attachment and Tool	15
Straight Attachments.....	15
Specialized Attachments.....	15
Activating the Motor	16
Disassembling the MR7 System.....	17
Depressurize the System	17
Disconnect Hoses	17
Discard the Lubricant/Diffuser Cartridge.....	17
Remove the Attachment and Tool.....	17
Cleaning & Sterilizing the MR7 System.....	18
MR7 Motor	18
Legend Attachments / Tubes.....	19

MR7 Pneumatic Control Unit / Regulator Hose /
Triton Adapter / Instrument Case..... 21

Transmissible Spongiform Encephalopathies (TSE) Return Policy..... 21

Troubleshooting..... 22

Refurbishing or Repairs 25

Preventative Maintenance..... 26

Limited Warranty..... 27

Appendix A—Specialized Attachments..... 28

 Angled Attachments..... 28

 Angled Double Lock Attachments 29

 Curved Bur Attachments..... 30

 Variable Exposure Attachments 31

 Footed Attachments 32

 Rotating Footed Attachments 33

 Contra-Angle Attachment (16-MF)..... 34

 Metal Cutting Attachments..... 35

 Telescoping Attachments..... 36

 Perforator Driver Attachments..... 37

 Jacobs® Chuck Attachments 38

 Bone Mill Attachment..... 39

General Information

Read and understand this manual before use of the MR7 System.

The Midas Rex® MR7 system is designed for use by medical professionals familiar with powered surgical instrumentation. The surgeon is responsible for learning the proper techniques in the use of this system, as inappropriate use may potentially be harmful. It is strongly recommended that the surgeon and dedicated operating room personnel are knowledgeable with the use of this equipment by being trained in Medtronic Midas Rex Hands-On Workshops or by one of the local authorized representatives.

The MR7 system consists of the following components:

- MR7 or MR7 Touch Motor
- MR7 Pneumatic Control Unit with Various Connectors
- MR7 Regulator Hose
- MR7 Lubricant/Diffuser Cartridge
- MR7 Triton® Adapter (optional)
- Legend® Attachments*
- Legend® Dissecting Tools*

*The MR7 system uses the same attachments and dissecting tools as the Legend® Pneumatic High-Speed System.

Indications for Use

The Medtronic Midas Rex MR7 System is a pneumatically operated surgical instrument system. The pneumatic motors provide power to operate removable rotating surgical cutting tools and their accessories intended for use in neurosurgery, including craniotomy and spinal surgery; as well as Ear, Nose and Throat (ENT), orthopedic and general surgical applications including maxillofacial, craniofacial and sternotomy surgeries.

Contraindications

None

Special Notices

The words warning, caution and note have special meanings in this manual, and should be carefully reviewed:

WARNING: A warning indicates that the personal safety of the patient or physician may be involved. Disregarding this information could result in injury to the patient or physician.

CAUTION: A caution indicates that there is a risk of damaging equipment.

NOTE: A note is intended to provide additional information, which may be useful, but is not essential to complete the procedure.

General Safety Precautions

WARNINGS:

- Do not use the Midas Rex MR7 System before proper cleaning and sterilization.
- Do not operate the Midas Rex MR7 System in the presence of Magnetic Resonance Imaging devices.
- Do not use damaged, faulty, or modified Midas Rex MR7 System components. Inspect the Midas Rex MR7 System for damage prior to each use:
 - Check the motor's exhaust hose for cracks or tears.
 - Visually inspect attachments and tools. Do not use bent or damaged tools.
 - Install attachment and dissecting tool, then briefly run motor.
 - * Check motor for overheating and leaking lubricant.
 - * Check attachment for overheating.
 - * Check dissecting tool for flail.
- Do not operate the Midas Rex MR7 System without eye protection.
- Motors and attachments which fail due to extended use may allow a component to detach and fall from the motor or attachment, and may cause patient injury.
- Heavy side loads and/or long operating periods may cause the device to overheat. If overheating occurs:
 - Never place an overheated motor on the patient or draping during surgery.
 - Discontinue use and rest the motor by using intermittently, or wrap the motor/attachment interface with a moist sterile towel.
 - If the motor is passed off, the receiver should grasp the motor by the proximal end close to the motor hose.
- To avoid injury to the patient or user, do not place the handpiece on the patient or in an unsecured location, when not in use.
- Midas Rex MR7 motors should only be operated when the attachment is in the  position.

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High Vacuum Steam: 270° F (132° C) for 5 minutes

Gravity Displacement: 270° F (132° C) for 15 minutes

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No Latex Policy

Legend and MR7 products, packaging materials, labels, package inserts, and similar items manufactured by and/or for Medtronic Powered Surgical Solutions (MPSS) do not contain latex.


System Components

Non-Disposable Components

MR7 Motor

The MR7 motor is a high-speed, high-torque motor used to dissect bone and biomaterials.

Figure 1: Motor Components

1. Collet
2. Motor Case
3. Swivel
4. Hose
5. Finger Control Lever
6. Safety Slide (MR7 Touch Only)

In addition to the components listed in Figure 1, each MR7 motor has a lubricant/diffuser housing at the end of the motor hose, as seen in Figure 2.

Figure 2: Lubricant/Diffuser Housing

WARNING: Use only Medtronic Midas Rex Legend or MR7 devices with an MR7 motor. Use of other devices may cause injury or damage equipment, and will void the manufacturer's warranty.

Pneumatic Control Unit

The pneumatic control unit (Figure 3) provides variable speed motor control controls through a foot pedal. It also allows the user to switch between finger and foot control of the motor (if applicable).

Regulator

The regulator (Figure 4) controls the delivery pressure of compressed gas to the pneumatic control unit. The pressure gauges monitor cylinder pressure (right gauge) and delivery pressure (left gauge).

Note: Outlet pressure gauge accurate to +/- 12 psi.

Instrument Case

The instrument case (Figure 5) is used to organize equipment.

Regulator Hose

Connects from the gas source to the pneumatic control unit to deliver compressed gas.

Figure 3: Pneumatic Control Unit


Figure 4: Regulator


Figure 5: Instrument Case

N2 DISS to Male Schrader Adapter

The N2 DISS to male Schrader adapter (Figure 6) allows for the regulator hose to be attached to a female Schrader in-house gas connection. A N2 DISS to female Schrader adapter is also available for connection of the regulator hose to a male Schrader in-house gas connection.

Figure 6: N2 DISS to Male Schrader Adapter

N2 DISS to Air DISS Adapter

The N2 DISS to air DISS adapter (Figure 7) allows for the regulator hose to be attached to an Air DISS in-house gas connection.


Figure 7: N2 DISS to Air DISS Adapter

N2 DISS to WF4 Adapter

The N2 DISS to WF4 adapter (Figure 8) allows for the regulator hose to be attached to a Midas Rex safety valve regulator previously used for Midas Rex Classic or Midas Rex III motors. The in-line oiler must be removed from the safety valve regulator.


Figure 8: N2 DISS to WF4 Adapter

Motor Wrench

The motor wrench (Figure 9) is used to align arrows on motor collet flats prior to installation of a Legend attachment.


Figure 9: Motor Wrench

Triton Adapter

The Triton adapter (Figure 10) allows the Triton handpiece to be driven by the MR7 pneumatic control unit. It functions much the same way as the Triton port on the Legend pneumatic control unit, except that it is connected between the control unit and the gas source, rather than being integrated into the control unit.

Figure 10: Triton Adapter

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NOTE: Angled and straight attachments with the same length, marking, and color band share the same dissecting tool. Curved and straight telescoping tubes with the same length, marking, and color band also share the same dissecting tool. Example: The 14-AM straight and 14-AM angled attachments are 14 cm long, marked 14-AM and have a green color band. All dissecting tools with the prefix 14 (14MH30) may be used in either the 14-AM straight or 14-AM angled attachment.

Be sure to match the color code and nomenclature on the Legend Dissecting Tool packaging with the color band and nomenclature on the Legend Attachment.

Disposable Components

WARNING: Use only Medtronic Midas Rex Legend or MR7 devices with an MR7 motor. Use of other devices may cause injury or damage equipment, and will void the manufacturer's warranty.

Lubricant/Diffuser Cartridge

The lubricant/diffuser cartridge (Figure 11) provides lubrication to the motor and filters oil from exhausted air.

Figure 11: Lubricant/Diffuser Cartridge

Telescoping Tubes

Telescoping tubes (Figure 12) provide support to the rotating dissecting tool. Telescoping tubes are disposable following multiple uses and should be discarded when heat or excessive vibration is noticed or insertion of tools becomes difficult.

Figure 12: Telescoping Tube


Cleaning Brushes

Cleaning brushes (Figure 13) are used to clean debris from lumen of attachments and telescoping tubes. Sized for an internal bore diameter of 3.2 mm, 2.4 mm or 1.2 mm in Legend Attachments and Telescoping Tubes.

NOTE: Cleaning brushes will not pass through angled, contra-angle, metal cutting, perforator, or Jacobs Chuck attachments, because they are not cannulated.

Figure 13: Cleaning Brushes


Legend® Dissecting Tools

Legend dissecting tools are sterile cutting tools, intended for cutting bone and biomaterials.

Dissecting Tool Nomenclature

Part numbers for Legend dissecting tools follow a standard naming convention, which is described in the diagram below. A basic part number consists of the associated attachment length, the tool head shape, and the tool head diameter. Part numbers may also include a variety of prefixes to identify specific attachment types, as well as a variety of suffixes to provide additional information about the dissecting tool. Tools that use a design taken from the Mednext line are designated by an additional “-MN” suffix.

Tool Number Prefixes (not all inclusive)

F...	For use with footed attachments
MC	For use with metal cutting attachments
T	For use with telescoping attachments

Tool Head Shapes (not all inclusive)

AC	Acorn	MH	Match Head
BA	Ball	OV	Oval
CY	Cylinder	RT	Reverse Taper
HM	Hole Maker	TA	Tapered
HS	Hole Saw	TD	Twist Drill

Tool Number Suffixes (note that more than one of the suffixes listed may be combined in a single part number)

L	Long	S	Spiral
D	Diamond	SH	Short
X	Extra	DC	Diamond Coarse
F	Fine	DX	Diamond Extra Coarse
C	Carbide	MN	Mednext Tool Design


WARNINGS:

- Dissecting tools are for single-use only. Do not attempt to sterilize them. The dissecting tools are packed sterile and are not intended for repeat use. To prevent contamination, use only once.
- Do not use an attachment and dissecting tool combination that results in tool flail or excessive vibration.
- Do not attempt to remove a tool while the motor is running.
- Do not attempt to remove a tool from an overheated motor or attachment.
- Do not use the device if the package is opened or damaged.

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

WARNING: Be sure to match the color code and nomenclature on the Legend Dissecting Tool packaging with the color band and nomenclature on the Legend Attachment. Failure to do so could result in injury to the patient or operating room staff.

Surgical Application	Commonly Used Attachments	Commonly Used Dissecting Tools	Suggested Motor(s)
Spine	8-B, 9-M, 14-AM, 15-A	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection, etc.	MR7, MR7 Touch
		Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, sinus dissection, etc.	
		Oval  Helical cutting flutes and curved design blend acorn and ball styles to vary dissection efficiency with approach angle. For decortication, laminotomy, entry hole, nerve decompression, osteophyte removal, etc.	
		Hole Maker/Saw  Matched sets of Hole Makers and Hole Saws are efficient and effective for interbody fusion.	
		Cylinder  Effective bone sculpting and planing. For graft shaping, debridement, corpectomy, decortication, interbody fusion, fusion takedown, etc.	
		Acorn  Curved design varies dissection efficiency with varied approach angles. For entry hole, laminotomy, bone shaping, debridement, corpectomy, decortication, fusion takedown, etc.	
Telescoping		Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection, etc.	MR7, MR7 Touch
Footed, Straight		Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, entry hole, suture hole, midface advancement, etc.	

Surgical Application	Commonly Used Attachments	Commonly Used Dissecting Tools	Suggested Motor(s)
Neurosurgical–Cranial	7-6ST, 8-B, 9-M, 10-9ST, 14-AM, 15-A	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection, etc. Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, sinus dissection, etc. Twist Drill  Helical design with stop produces a hole with a precise depth. Ideal for plating. Acorn  Curved design varies dissection efficiency with varied approach angles. For entry hole, laminotomy, bone shaping, debridement, corpectomy, decortication, fusion takedown, etc.	MR7
	Telescoping	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection, etc.	
	Footed	Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, entry hole, suture hole, midface advancement, etc.	
General Surgery and Plastic Surgery (Craniofacial/Maxillofacial/Sternotomy)	7-6ST, 8-B, 9-M, 10-9ST, 14-AM	Match Head  Elongated spherical design allows controlled, delicate dissection. For entry hole, nerve decompression, osteophyte removal, sinus dissection, etc. Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, sinus dissection, etc. Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, entry hole, suture hole, midface advancement, etc. Twist Drill  Helical design with stop produces a hole with a precise depth. Ideal for plating.	MR7, MR7 Touch
Ear, Nose, and Throat (Otology, Neurotology)	7-6ST, 10-9ST	Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, sinus dissection, etc.	MR7

Surgical Application	Commonly Used Attachments	Commonly Used Dissecting Tools	Suggested Motor(s)
Orthopaedics	8-B, 9-M, 14-AM, 21-TU, 26-R, Footed, Telescoping	Ball  Helical cutting flutes dissect bone or cement effectively from a wide variety of approach angles. For debridement, decortication, sinus dissection, etc.	MR7, MR7 Touch
		Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, entry hole, suture hole, midface advancement, etc.	
		Acorn  Curved design varies dissection efficiency with varied approach angles. For entry hole, laminotomy, bone shaping, debridement, corpectomy, decortication, fusion takedown, etc.	
		Cylinder  Effective bone sculpting and planing. For graft shaping, debridement, corpectomy, decortication, interbody fusion, fusion takedown, etc.	
Biometals/ Bioceramics/ Biomaterials	Footed	Tapered  Slender design for precise dissection with minimal bone loss. For transection, osteotomy, graft harvesting, bone shaping, entry hole, suture hole, midface advancement, etc.	MR7, MR7 Touch
	MC	Metal Cutter  Cutting flutes or diamond wheel design remove metals, ceramics and other biomaterials effectively from a variety of approach angles. For cutting rods, pins, plates, implants, screws, etc.	

Setting up the Operating Room

Power Source Requirements


Required Operating (Dynamic) Pressure	Nominal Operating (Dynamic) Pressure	Approximate Flow Rate Required	Gas Type
80–120 psi	100 psi	12 cubic feet/min.	Nitrogen or Dry-Filtered Compressed Air
5.5–8.3 bar	6.9 bar	340 liters/min.	

CAUTION: Do not run the motor at an operating pressure above or below the required operating pressure range. Operating pressure below 80 psi (5.5 bar) may not provide proper lubrication to the motor. Operating pressure above 120 psi (8.3 bar) may damage or reduce the life of the motor.

Pneumatic Connections

Figure 14: Gas Connection Options

- 1. Regulator Hose (N2 DISS)
- 2. N2 DISS to Air DISS Adapter Gas Source
- 3. N2 DISS to Male Schrader Adapter Gas Source
- 4. Gas Source (N2 DISS)
- 5. Regulator
- 6. DISS/WF4 Adapter
- 7. Regulator

CAUTION: If you are using the Midas Rex Safety Valve Regulator instead of the Legend Regulator, you must replace the in-line oiler with the DISS/WF4 adapter before use.

Setting up the MR7 System

Installing the Oiler Cartridge

WARNING: Do not use the MR7 system with the Midas Rex in-line oiler. The MR7 motor is sufficiently lubricated by the lubricant/diffuser on the motor hose, and will be over-lubricated if the Midas Rex in-line oiler is used.

1. Set the non-running (static) pressure to 80–120 psi (5.5–8.3 bar) at the gas source. Operating (dynamic) pressure may be adjusted later.
2. Hold the lubricant/diffuser cartridge perpendicular to the housing (Figure 15), and press the cartridge's circular fitting onto the housing's circular receptacle (Figure 16), breaking the foil seal.
3. Rotate the cartridge down until it clicks into place.
4. Verify that the  symbol on the cartridge is lined up with the notch on the housing (Figure 17).

WARNINGS:

- Failure to properly secure the lubricant/diffuser cartridge may cause injury to operator and/or operating room staff.
- Do not attempt to remove the lubricant/diffuser cartridge while the system is pressurized.

CAUTIONS:

- Do not use an MR7 motor without a lubricant/diffuser installed.
- Do not use a lubricant/diffuser cartridge for more than one hour of drill time.
- Do not re-use a lubricant/diffuser cartridge. It is a single-use product.
- Do not attempt to refill a used lubricant/diffuser cartridge.
- Do not use a lubricant/diffuser cartridge if it appears to be damaged, or if the inner foil seal is punctured.

Figure 15: Aligning the Lubricant/Diffuser Cartridge with the Housing



Figure 16: Pressing the Cartridge onto the Housing

Figure 17: Correctly Installed Lubricant/Diffuser Cartridge

Connecting the Motor

Connect the motor hose to the motor port on the top of the pneumatic control unit, by swinging the port cover to the side and pressing the end of the hose into the port (Figure 18).

WARNING: Do not pinch, kink, obstruct, cut, tear, or step on the motor/exhaust hose. This may cause the hose to burst, potentially injuring the patient or user.

NOTES:

- If using the MR7 Touch motor, slide the control slide on the pneumatic control unit to the  position (Figure 19). This will automatically depress and lock the foot pedal. The control will not lock into the  position unless the motor hose is connected into the motor port. When the motor hose is removed from the motor port, the foot pedal will return to normal position.
- If using the Triton Power Surgical Instrumentation System in conjunction with the MR7 motor, use the optional Triton adapter to connect the Triton hose. Refer to the documentation accompanying the adapter for connection instructions.
- The motor's exhaust hose may have an oily film on the external surface from pressure and/or temperature differentials following sterilization. Wipe the exhaust hose with a sterile cloth prior to use. If motor continues to have oil on the exhaust hose, return the motor to MPSS for refurbishing.

WARNING: To avoid injury to the patient or user, do not use the pneumatic control unit to operate systems other than the MR7, Legend, and Triton systems.

Prior to installation of a Legend attachment and dissecting tool, ensure that the arrows on the motor collet flats are aligned (Figure 20). If the arrows are not aligned, use the motor wrench to turn the collet flat closest to the motor case until its arrow is aligned with the arrow on the other collet flat.

WARNING: To avoid injury when using the MR7 Touch motor, ensure the safety slide is in the "O" position before installing the attachment and tool.

Figure 18: Connecting the Motor Hose to the Motor Port

Figure 19: Finger/Foot Control Slide



Figure 20: Aligning the collet flats

Installing an Attachment and Tool

WARNING: Dissecting tool flutes are sharp and may perforate surgical gloves. Tools may be grasped with a hemostat to aid in installation and removal.

Straight Attachments

Installation:

1. Slide the attachment over the motor collet, aligning the triangular markers (Figure 21). You will feel and hear the attachment click into place when it is fully seated.
2. Insert the dissecting tool into the attachment with a slight rotational motion (Figure 22). You will feel and hear the tool click into place when it is fully seated in the attachment.
3. Turn the attachment to the  position on the motor case (Figure 23). Gently pull on the shaft of the dissecting tool to verify proper installation.

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Specialized Attachments

See *Appendix A—Specialized Attachments* for installation and removal instructions for other attachments.

Figure 21: Sliding the Attachment over the Motor Collet


Figure 22: Inserting the Dissecting Tool into the Attachment



Figure 23: Attachment in the Locked Position

Activating the Motor

NOTE: In order to activate the MR7 Touch motor, the safety slide on the finger control switch must be in the | position, and the control slide on the foot control must be in the  position. The control slide will not lock in the  position unless a motor is connected to the motor port.

1. Activate the motor by pressing on the foot control pedal (Figure 24), or by pressing on the finger control lever (MR7 Touch motor only).
2. Adjust operating pressure as needed at the compressed gas source until supply pressure gauge on pneumatic control unit reads within a range of 80–120 psi (5.5–8.3 bar) as required. Operating pressure (with motor running) will decrease slightly from the non-running (static) pressure setting when the motor is activated.

Figure 24: Foot Control Pedal


Figure 25: MR7 Finger Control Switch

WARNINGS:

- Do not use excessive force to pry or push bone with the attachment or tool during dissection. This could cause the tool to break and cause injury to the patient or operating room staff.
- Use adequate irrigation during dissection, to prevent thermal necrosis.
- MR7 motors should only be operated when the attachment is in the  position.

NOTE: To decrease pressure, turn down the in-house compressed gas source or loosen the pressure handle on the regulator. Push down on the pressure relief at the pneumatic control unit to exhaust excess pressure in the hoses. Then re-adjust pressure as needed.

Disassembling the MR7 System

Depressurize the System


1. Turn off the compressed gas at the source.
2. Press the pressure relief button on the pneumatic control unit, to release remaining gas.

WARNING: Do not disassemble equipment before the gas is released from control unit.

Disconnect Hoses

Release the motor hose from the control unit, by holding the hose firmly and pressing the motor port connection button.

Figure 26: Motor Port Connection and Pressure Relief Buttons

Discard the Lubricant/Diffuser Cartridge

Remove the lubricant/diffuser cartridge from the housing and discard it.

CAUTION: Do not re-use a lubricant/diffuser cartridge. It is a single-use product.

Remove the Attachment and Tool

1. Follow the applicable removal instructions in the *Installing an Attachment and Tool* section of this manual to remove the attachment and tool from the motor.

WARNING: Dissecting tool flutes are sharp and may perforate surgical gloves. Tools may be grasped with a hemostat to aid in installation and removal.

2. Discard used dissecting tools in an appropriate container.

Cleaning & Sterilizing the MR7 System

MR7 Motor

Warnings and Precautions	Do not soak/submerge MR7 devices. Do not use ultrasound to clean MR7 devices. Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. Allow an adequate cooling period after steam sterilization. Use ONLY nylon cleaning brushes. Non-nylon cleaning brushes leave residue that may prevent the tool from being secured properly in the handpiece.				
Limitations	Verify functionality prior to re-use.				
INSTRUCTIONS					
Point of Use	No particular requirements. Follow hospital procedures.				
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination	No particular requirements. Follow hospital procedures.				
Cleaning: Automated (Do NOT use ultrasonic washer)	Review the washer-disinfector warning above, before using this cleaning method. Remove devices from instrument trays before placing into washer baskets. Orient devices following recommendations of the washer/disinfector manufacturers. Verify that devices are visually clean after automated cleaning. Recommended Washer Cycle Pre-Wash: Cold tap water, 2 min. Wash: 66oC, 5 min. using a neutral enzymatic detergent, pH 6.0–8.0 Rinse: Hot tap water, 1 min.				
Cleaning: Manual	Wipe all external surfaces of the motor and hose, and wipe inner surface of oiler housing with a cloth dampened with a neutral enzymatic detergent, pH 6.0–8.0. Brush motor case and collet with a nylon brush dampened with a neutral enzymatic detergent. If using the MR7 Touch motor, be sure to brush under the finger control lever. Rinse motor thoroughly under running water, collet end pointed down. Dry collet and motor with lint free towel. Verify that devices are visually clean after manual cleaning.				
Disinfection	No particular requirements. Follow hospital procedures.				
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap				
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization:				
	Cycle:	Gravity	Pre-vac	Pre-vac (FR/WHO)	Pre-vac (UK)
	Temperature:	132°C	132°C	134°C	134°C
	Time:	25 min.	4 min.	18 min.	3 min.
	Drying:	15 minutes	15 minutes	20 minutes	10 minutes
	*Items contaminated with Transmissible Spongiform Encephalopathies (TSE) agents may be decontaminated using steam autoclave at a temperature of 134–137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles as referenced in NHS Estates HTM 2010 parts 4 & 6: Appendix 2, Items Contaminated With TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies. Medtronic recommends that all products used on a patient suspected or confirmed with a TSE diagnosis be incinerated.				
	STERRAD Sterilization: Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.				
	100% EtO Sterilization Parameters: Preconditioning: 51–59°C, 70 ±5% relative humidity, 30 min. Temperature: 51–59°C Relative humidity: 70 ±5% Ethylene oxide concentration: 725 ± 25mg/L Gas exposure time (full-cycle): 4 hours Aeration: 18 hours at 51–59°C				
Steris: Do not use liquid peracetic acid sterilization due to immersion procedure.					

Maintenance, Inspection and Testing	Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. After cleaning and sterilization, verify functionality prior to re-use.
Storage	Store with other sterile devices.

NOTE: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

Legend Attachments / Tubes

Warnings and Precautions	Do not soak/submerge Legend devices. Do not use ultrasound to clean Legend devices. Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. Allow an adequate cooling period after steam sterilization. Use ONLY nylon cleaning brushes. Non-nylon cleaning brushes leave residue that may prevent the tool from being secured properly in the handpiece.
Limitations	Verify functionality prior to re-use.

INSTRUCTIONS	
Point of Use	No particular requirements.
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.
Preparation for Decontamination	No particular requirements. Follow hospital procedures.
Cleaning: Automated (Do NOT use ultrasonic washer)	Review the washer-disinfector warning above, before using this cleaning method. Manually rinse attachments/tubes under tap water, until no visible soil is noticed, before placing them into the automatic washer. Remove devices from instrument trays before placing into washer baskets. Orient devices following recommendations of the washer/disinfector manufacturers. Recommended Washer Cycle Pre-Wash: Cold tap water, 2 min. Wash: 66°C, 5 min. using a neutral enzymatic detergent, pH 6.0–8.0 Rinse: Hot tap water, 1 min.
Cleaning: Manual	Wipe all attachments and telescoping tubes with a cloth, dampened with a surgical instrument cleaning solution. Immerse the head of Contra-Angle attachments in surgical instrument cleaning solution and run the motor for 1 minute. Other attachments and tubes may be mechanically agitated in cleaning solution, but not soaked or immersed. A nylon brush dampened with a surgical instrument cleaning solution may be used to clean the external surfaces and internal connecting surfaces of the attachments and tubes. Straight attachments, footed attachments and telescoping straight tubes have special cleaning brushes sized to the attachment's or telescoping tube's internal diameter. Push the brush wet with surgical instrument cleaning solution through the attachment or telescoping tube from rear to front to loosen and remove debris trapped inside. Move any moveable parts back and forth to allow solution to thoroughly clean attachment, e.g., sleeve on footed attachment, perforator attachment. Rinse thoroughly with tap water. Thoroughly dry attachments. An air gun may be used to blow moisture out from rear to front of attachment. Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: Holding the can approximately 10–15 cm (3–6 in.) away from the attachment, spray all components that move, rotate, or slide with three quick squirts. Articulate movable components to ensure proper lubrication. Remove excess lubricant with a clean cloth.
Disinfection	No particular requirements Follow hospital procedures.
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap.

Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization:				
	Cycle:	Gravity	Pre-vac	Pre-vac (FR/WHO)*	Pre-vac (UK)*
	Temperature:	132°C	132°C	134°C	134°C
	Time:	25 min.	4 min.	18 min.	3 min.
	Drying:	15 minutes	15 minutes	20 minutes	10 minutes
	*Items contaminated with TSE agents may be decontaminated using steam autoclave at a temperature of 134–137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles as referenced in NHS Estates HTM 2010 parts 4 & 6: Appendix 2, Items Contaminated With TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies. Medtronic recommends that all products used on a patient suspected or confirmed with a TSE diagnosis be incinerated.				
	STERRAD Sterilization: Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.				
	100% EtO Sterilization Parameters: Preconditioning: 51–59°C, 70 ±5% relative humidity, 30 min. Temperature: 51–59°C Relative Humidity: 70 ±5% Ethylene oxide concentration: 725 ± 25mg/L Gas exposure time (full-cycle): 4 hours Aeration: 18 hours at 51–59°C				
	Steris: Do not use liquid peracetic acid sterilization due to immersion procedure.				
Maintenance, Inspection and Testing	Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. Verify functionality prior to re-use.				
Storage	Store with other sterile devices.				

NOTE: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

MR7 Pneumatic Control Unit / Regulator Hose / Triton Adapter / Instrument Case

Warnings and Precautions	Do not soak/submerge MR7 devices. Do not use ultrasound to clean MR7 devices. Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. Do not sterilize the MR7 pneumatic control unit, regulator hose, or Triton adapter.
Limitations	Verify functionality prior to re-use.

INSTRUCTIONS	
Point of Use	No particular requirements. Follow hospital procedures.
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.
Preparation for Decontamination	No particular requirements. Follow hospital procedures.
Cleaning: Automated (Do NOT use ultrasonic washer)	Not validated.
Cleaning: Manual	Wipe the pneumatic control unit, regulator hose, Triton adapter, and instrument case with a cloth dampened with surgical instrument cleaning solution after each use.
Disinfection	No particular requirements. Follow hospital procedures.
Packaging	No particular requirements. Follow hospital procedures.
Sterilization	Do not sterilize pneumatic control unit, regulator hose, or Triton adapter.
Maintenance, Inspection and Testing	Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. Verify functionality prior to re-use.
Storage	Do not store with sterile devices.

NOTE: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

Transmissible Spongiform Encephalopathies (TSE) Return Policy

Medtronic Powered Surgical Solutions (MPSS) will not authorize or accept the return of MPSS products that directly contact patients or is contaminated with a patient's body fluids who is suspected or confirmed with a Transmissible Spongiform Encephalopathies / Creutzfeldt-Jakob Disease (TSE/CJD) diagnosis. Furthermore, MPSS recommends that all MPSS products used on a patient confirmed with a TSE diagnosis be incinerated. Contact your Sales Representative for replacement of product incinerated under this policy or for temporary equipment while original equipment is quarantined. Contact MPSS Regulatory Affairs Department for additional information regarding TSE contamination. MPSS dissecting tools used on a patient suspected of a TSE/CJD diagnosis must be incinerated.

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Troubleshooting

NOTE: All Legend and MR7 devices returned for servicing or refurbishing should be properly cleaned and sterilized prior to shipping.

Motor not Running or Low on Power:

Possible Cause	Solution
Hoses not properly connected.	Make sure all connections are secure.
Operating pressure inadequate.	Check gas supply pressure gauge. Increase pressure according to compressed gas requirements, if necessary.
Attachment not properly installed and locked onto the motor.	Remove and re-install attachment and tool to ensure proper installation and locking of attachment onto motor.
Foot pedal on pneumatic control unit not functioning properly.	Check for obstructions under the foot pedal. If foot pedal continues to fail, return the pneumatic control unit to MPSS to be refurbished.
Motor stalls.	Manually spin the dissecting tool, then activate the motor. If the motor continues to stall, return it to MPSS to be refurbished.

Motor Continues to Run:

Possible Cause	Solution
Pneumatic control unit is not functioning properly.	Depressurize the system and return the pneumatic control unit to MPSS to be refurbished.
Finger control is not functioning properly.	Return motor to MPSS to be refurbished.
Pneumatic control unit is locked in the finger control position.	Move the finger control lever to the foot control position.

System Makes an Abnormal Noise:

Possible Cause	Solution
Inadequate lubrication.	Check for proper installation of the lubricant/diffuser cartridge. If the problem persists, replace the cartridge. If replacing the cartridge doesn't fix the problem, return the motor to MPSS for refurbishing.
Motor's exhaust hose is damaged, or internal pressure hose is detached.	Depressurize the system and return the motor to MPSS to be refurbished.
Worn bearings.	Switch attachments to determine whether the bearings are failing in the motor or in the attachment. Return the failing component to MPSS to be refurbished.
Attachment not properly installed and locked onto the motor.	Remove and reinstall attachment and tool to ensure proper installation and locking of the attachment onto the motor.
Safety relief valve has been activated by high air pressure.	Ensure that air operating/dynamic air pressure is no higher than 120 psi.

Motor is Too Hot to Touch/Hold:

Possible Cause	Solution
Inadequate cool down period following sterilization.	Motor must be allowed to cool down following steam sterilization.
Inadequate lubrication.	Check for proper installation of the lubricant/diffuser cartridge.
Attachment transferring heat to the motor.	Switch attachments to determine whether the heat is being generated by the motor or the attachment. Return the failing component to MPSS for refurbishing.
Heavy side loading during dissection.	Discontinue use and rest the motor by using it intermittently or wrap the motor with a moist sterile towel. If overheating continues, return the motor to MPSS for refurbishing.
Inadequate irrigation.	Ensure adequate irrigation to surgical site during bone dissection.

Attachment Will not Properly Seat on the Motor:

Possible Cause	Solution
Motor collet flats are not aligned.	Use the Legend motor wrench to rotate the flat closest to the motor case until its marker is aligned with the marker on the flat farthest away from the motor case.

Tool is Difficult to Remove from Attachment:

Possible Cause	Solution
Aging of attachment.	Return to MPSS to be refurbished, or purchase new equipment.
Improper cleaning.	
Use of reprocessed tools.	
Use of an unauthorized refurbisher.	

16-MF Contra-Angle Attachment is Overheating:

Possible Cause	Solution
The Contra-Angle attachment operates by a set of internal gears to engage the drive shaft. It is normal for some heat to be generated approximately 2cm from the distal end of the attachment and at the right of the angle head.	Verify pressure setting of 80 psi (5.5 bar). If heat continues or is excessive, return the attachment to MPSS to be refurbished or purchase new equipment.

Perforator is Running too Slow:

Possible Cause	Solution
Pressure set incorrectly.	Check the pressure setting at the foot control.

Dissecting Tool Flails:

Possible Cause	Solution
A non-Legend dissecting tool is being used.	Replace with a Legend dissecting tool.
Worn attachment or tube bearings.	Try another attachment or tube to isolate the location of the problem. If the attachment is failing, return it to MPSS. If the tube is failing, dispose of it and use a new tube.
Attachment/tube and tool are not compatible.	Match color code on the dissecting tool packaging to the color code on the attachment/tube.
Motor is damaged.	Return motor to MPSS to be refurbished.
Dissecting tool's size and geometry may contribute to flailing at certain speeds.	Adjust the speed by changing the pressure setting or foot/finger control. Do not use if flailing persists. Change dissecting tools.

Dissecting Tool Vibrates Excessively:

Possible Cause	Solution
Dissecting tool's size and geometry may create excessive vibration at certain speeds.	Adjust the speed by changing the pressure setting or foot/finger control. Change dissecting tools.

Dissecting Tool Will not Seat Properly in the Motor or Attachment Collet:

Possible Cause	Solution
Debris in collet of attachment or motor.	Clean the attachment or motor thoroughly according to the instructions in this manual. If cleaning does not correct the problem, return the attachment or motor to MPSS to be refurbished.
A non-Legend dissecting tool is being used.	Replace with a Legend dissecting tool.

Smoke is Generated by the Attachment or Motor:

Possible Cause	Solution
Attachment is not in the locked position.	Make sure the attachment is in the locked position.

Refurbishing or Repairs

When the MR7 System requires servicing or refurbishing, contact Medtronic Powered Surgical Solutions Repair Services for a return authorization and instructions for returning the equipment. Medtronic Powered Surgical Solutions provides quality assured service by factory-trained personnel who will utilize genuine Midas Rex Legend parts as required. All items being returned for servicing or refurbishing should be properly cleaned and sterilized prior to shipping.

Contact:

Medtronic Powered Surgical Solutions Repair Services:
(800) 335-9557 or (817) 788-6440
mmrcustomerservice@medtronic.com

Peak performance, reliability and maximum service life from your MR7 System may be assured by using only those Midas Rex Legend products for your MR7 System that are manufactured by and sold through Medtronic Powered Surgical Solutions, Fort Worth, Texas. While Medtronic Powered Surgical Solutions guarantees complete compatibility among its products within a specific product line, the dissecting tools are designed for single-use only, and Medtronic disclaims any responsibility when reprocessed dissecting tools are used. If you would like more information about the patient and product risks associated with reprocessed tools, please contact the number or e-mail address listed above.

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Preventative Maintenance

The *Midas Rex MR7 System Preventive Maintenance Manual* has been developed to assist you in getting the greatest ownership value from your MR7 System, while helping to maximize its performance, safety and reliability. The scheduled preventive maintenance/service program is in addition to the required routine cleaning after each use. Please refer to the preventive maintenance manual for the specific steps necessary to maintain the MR7 System.

Limited Warranty

- A.** This Limited Warranty provides the following assurance to the purchaser of a Medtronic Midas Rex® MR7 Pneumatic High Speed System. This Limited Warranty is extended only to the buyer purchasing the MR7 System directly from Medtronic or from its affiliate or its authorized distributor or representative. The Midas Rex® MR7 Pneumatic High Speed System includes the motor, foot control, instrumentation cases and trays (hereafter referred to as System Components), straight and angled motor attachments (hereinafter referred to as "Attachments"), telescoping tubes (hereinafter referred to as Semi-reusable Components) and dissecting tools and other accessories not listed above and jointly referred to as MR7 Pneumatic High Speed System, unless specifically noted.
- (1) Should a System Component fail to function to Medtronic's published specifications during the term of this Limited Warranty (one year from the date of sale of a new System Component or 90 days from the date of sale of a refurbished or used System Component), Medtronic will either repair or replace the Motor Component or any portion thereof.
 - (2) Should an Attachment fail to function to Medtronic's published specifications during the term of this Limited Warranty (90 days from the date of sale of a new Attachment), Medtronic will either repair or replace the Attachment or any portion thereof.
 - (3) Should a Semi-reusable Component fail to function to Medtronic's published specifications during the term of this Limited Warranty (30 days from the date of sale of a new Semi-reusable Component), Medtronic will replace the Semi-reusable Component or any portion thereof.
 - (4) Should a Single Use Component fail to function to Medtronic's published specifications prior to its "use by" date Medtronic will replace the Single Use Component.
- B.** To qualify for this Limited Warranty, these conditions must be met:
- (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
 - (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
 - (3) Medtronic must be notified in writing within thirty (30) days following discovery of a defect.
 - (4) The Product must be returned to Medtronic within thirty (30) days of Medtronic receiving notice as provided for in (3) above.
 - (5) Upon examination of the Product by Medtronic, Medtronic shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services, if applicable, have been performed on the Product
- C.** This Limited Warranty is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the MR7 System, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
- D.** The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this Limited Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

Appendix A—Specialized Attachments

Angled Attachments

Installation:

- 1. Slide the attachment over the motor collet, aligning the triangular markers (Figure 27). You will feel and hear the attachment click into place when it is fully seated.
- 2. Turn the attachment to the  position on the motor case (Figure 28).
- 3. With the “tool lock” in the  position, insert the dissecting tool into the attachment with a slight rotational motion (Figure 29). You will feel and hear the tool click into place when it is fully seated.
- 4. Turn the tool lock to the  position (Figure 30). Gently pull on the shaft of the dissecting tool to verify proper installation.

NOTE: A dissecting tool may be installed in the attachment before the angled attachment is placed on the motor.

CAUTION: Hold the handpiece assembly by the attachment, so that the attachment does not inadvertently loosen from the handpiece.

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Figure 27: Sliding the Attachment over the Motor Collet


Figure 28: Attachment in the Locked Position


Figure 29: Inserting the Dissecting Tool into the Attachment



Figure 30: Tool Lock in the Locked Position

Angled Double Lock Attachments

Installation:

1. Slide the attachment over the motor collet, aligning the triangular markers.
2. Pull the attachment towards the motor and turn the attachment to the  position.
3. With the “tool lock” in the  position, insert the dissecting tool into the attachment with a slight rotational motion. You will feel the tool click into place when it is fully seated.
4. Turn the tool lock to the  position. Gently pull on the shaft of the dissecting tool to verify proper installation.
5. On Variable Exposure Attachments, use the TUBE adjustment ring to adjust the exposure of the dissecting tool. With the tool pointing away from you, turn the ring to the right to increase the length of the tube, thereby decreasing the exposure of the tool. Turn the ring to the left to decrease the length of the tube, thereby increasing the exposure of the tool.

Removal:

1. To remove the attachment, hold the motor in the palm of your hand, and push the sleeve on the attachment distally while turning the attachment to the  position.
2. Release the sleeve and remove the attachment.

Figure 31: Installing the attachment.

Figure 32: Inserting the Dissecting Tool into the Attachment

Figure 33: Removing the Attachment

Curved Bur Attachments

WARNINGS:

- Use adequate irrigation and keep the cooling sleeve soaked during dissection. Inadequate irrigation may cause thermal necrosis.
- Do not modify the bur. Bending or prying may break the bur, causing harm to the patient or operating room staff.
- Excessive pressure applied to the bur may cause bur damage. If this occurs, use extreme care to ensure that all fragments of the bur are removed from the patient.
- Disposable devices are for single-use only. Do not attempt to sterilize disposable devices. The disposables are packed sterile and are not intended for repeat use. To prevent contamination, use only once.
- Test for bur wobble (eccentricity) at the desired speed prior to use. Select a new bur or reduce speed if wobble is observed prior to use or during the procedure. Bur wobble may cause patient injury.

CAUTIONS:

- Curved burs are not designed for variable tool exposure. Do not attempt to adjust tool exposure, as this may damage the device.
- When operating or testing the motor, ensure the bur is properly inserted and locked into AT10 attachment.

Installation

1. Turn the tube locking ring and the tool locking ring on the AT10 (V01) Telescoping Attachment into the  position.
2. Slide the curved bur into the AT10 (V01), until the hub is fully seated.
3. Gently press on the head of the dissecting tool to lock it into place. You will hear a click when the curved bur is fully seated.
4. Turn both locking rings to the  position on the AT10 (V01).
5. Gently pull on the curved bur to verify proper installation.
6. Refer to your motor's *Instruction Manual* for information on setting up the AT10 (V01) attachment and motor.
7. Briefly run the motor with the curved bur installed, checking for bur wobble or excessive vibration. If either occurs, perform the corrective actions described in the *Troubleshooting* section of these instructions.

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Removal

1. Turn both AT10 (V01) locking rings to the  position and pull the bur out of the attachment.
2. Discard the curved bur, according to hospital procedures.

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Variable Exposure Attachments

CAUTION: The Legend Variable Exposure attachments can be distinguished from standard attachments by the dual color bands on the attachment. Match the color band on the attachment to the color code on the dissecting tool packaging.

WARNINGS:

- Surgeons should familiarize themselves with the performance of dissecting tools before use, and should explore the effect of various levels of tool exposure on dissection stability. If the tool exhibits excessive chatter, vibration, or movement, decrease the tool exposure.
- Dissecting tool size and geometry may contribute to excessive vibration at certain speeds. Increase or decrease speed by adjusting the foot/finger control, or by changing the operating pressure or console speed setting. If necessary, use a different dissecting tool.

Installation:

1. Refer to the appropriate section of this manual on installing fixed or angled attachments.
2. After installation, use the TUBE adjustment ring to adjust the exposure of the dissecting tool (Figure 34). With the tool pointing away from you, turn the ring to the right to increase the length of the tube, thereby decreasing the exposure of the tool. Turn the ring to the left to decrease the length of the tube, thereby increasing the exposure of the tool.

WARNING: Do not use the Variable Exposure Attachment if the TUBE adjustment ring spins freely or fails to click into place with each adjustment, as the exposure may change without warning.


CAUTION: Make sure that the tool lock (angled attachments only) and the attachment lock are still in the  position after each adjustment of the tool exposure. Attempting to increase the exposure too far may result in the attachment becoming unlocked. Accidentally turning the tool lock may result in reduced speed and/or overheating of the attachment.

WARNING: Do not use the end of the tube as a depth gauge or depth stop.

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Footed Attachments

Installation:

- 1. Insert the dissecting tool into the motor collet with a slight rotational motion (Figure 35). You will feel and hear the tool click into place when it is fully seated.
- 2. Slide the footed attachment over the dissecting tool, onto the motor, aligning the triangular markers (Figure 36).
- 3. Pull the footed attachment towards the motor and turn the attachment to the  position (Figure 37).

Removal:

- 1. To remove the Legend footed attachment, hold the motor in the palm of your hand, and push the sleeve on the attachment distally while turning the attachment to the  position (Figure 38).

WARNING: Remove Legend Footed Attachments cautiously and slowly per instructions to avoid injury to the operator.

- 2. Release the sleeve and remove the attachment.
- 3. Pull the dissecting tool out of the motor collet and discard the tool.

Figure 35: Inserting the Tool into the Collet


Figure 36: Sliding the Attachment onto the Motor

Figure 37: Attachment in the Locked Position


Figure 38: Removing the Attachment

Rotating Footed Attachments

NOTE: Rotating and fixed footed attachments with the same length, marking, and color band share the same dissecting tools.

Installation:

1. Insert the dissecting tool into the motor collet with a slight rotational motion (Figure 39). You will feel and hear the tool click into place when it is fully seated.
2. Slide the attachment over the dissecting tool and onto the motor, aligning the triangular markers (Figure 40). You will feel and hear the attachment click into place when it is fully seated.
3. Turn the attachment to the  position (Figure 41).

NOTE: The footed end of the attachment now has 360° of unrestricted rotation.

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Figure 39: Inserting the Tool into the Collet

Figure 40: Sliding the Attachment onto the Motor



Figure 41: Attachment in the Locked Position

Contra-Angle Attachment (16-MF)

CAUTION: Do not run the 16-MF attachment with operating pressure above 80 psi (5.5 bar). This may cause over-heating and damage to internal gears of attachment.

Installation:

1. Decrease the pressure at the compressed gas source to 80 psi (5.5 bar).
2. Adjust the pressure as needed by lowering it at the gas source, then push down the pressure relief on the pneumatic control unit to exhaust the excess pressure in the hoses.
3. Slide the attachment over the motor collet, aligning the triangular markers (Figure 42). You will feel and hear the attachment click into place when it is fully seated.
4. Turn the attachment to the  position (Figure 43).
5. Turn the attachment head lever to the open position (Figure 44).
6. Insert the dissecting tool and return the lever to the closed position. Gently pull on the dissecting tool shaft to verify proper installation.

Removal:

1. Removal is the reverse of installation.
2. Discard the dissecting tool after removing it from the attachment.

NOTE: A dissecting tool may be installed and locked in the attachment before the Contra-Angle attachment is installed on the motor.

Figure 42: Sliding the Attachment over the Motor Collet

Figure 43: Attachment in the Locked Position


Figure 44: Opening the Attachment Head Lever

Metal Cutting Attachments

WARNING: For metal transection, observe the following safety guidelines:

- Wear eye protection.
- Irrigate well to cool the cutting surfaces.
- Protect the wound site from metal debris.
- Use a clamp or grasping device to control loose fragments during transection of any metal component.

The Metal Cutting Attachment uses the tungsten carbide or diamond wheel dissecting tools.

Installation:

1. Slide the attachment over the motor collet aligning the triangular markers (Figure 45). You will feel and hear the attachment click into place when it is fully seated.
2. Turn the attachment to the  position.
3. Turn the tool lock towards the  icon, and insert the dissecting tool into the attachment with a slight rotational motion (Figure 46). You will feel and hear the tool click into place when it is fully seated.
4. Turn the tool lock to the  position (Figure 47). Gently pull on the shaft of the dissecting tool to verify proper installation.

NOTE: A dissecting tool may be installed and locked in the attachment before the metal cutting attachment is installed on the motor.

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WARNING: Do not use metal cutting dissecting tools on bone.

Figure 45: Sliding the Attachment over the Motor Collet

Figure 46: Inserting the Tool into the Attachment

Figure 47: Turning the Tool Lock to the Locked Position

Telescoping Attachments

Installation:

1. Slide the attachment over the motor collet aligning the triangular markers (Figure 48). You will feel and hear the attachment click into place when it is fully seated.

NOTE: A dissecting tool may be installed and locked in the attachment before the telescoping attachment is installed on the motor.

2. Turn the attachment to the  position.
3. Turn the TUBE locking ring towards the  icon.
4. Insert the base of the telescoping tube into the attachment (Figure 49).
5. To lock the tube in place, turn the TUBE locking ring towards the  icon. Do not over tighten.
6. To insert the tool, make sure that the TOOL locking ring is in the  position, and insert the dissecting tool into the top of the tube (Figure 50). You will feel and hear the tool click into place when it is fully seated.
7. Turn the TOOL locking ring to the  position. Gently pull on the shaft of the dissecting tool to verify proper installation.
8. If the tube position needs to be changed, turn the TUBE locking ring towards the  icon, re-position the tube, then turn the TUBE locking ring towards the  icon.

Removal:

1. To remove the attachment, unlock both locking rings, and pull the telescoping tube and tool out of the attachment.
2. Turn the attachment to the  position and remove it from the motor.

NOTE: Telescoping tubes are disposable following multiple uses, and should be discarded when heat or excessive vibration is noticed, or when insertion of the tool becomes difficult.

Figure 48: Sliding the Attachment over the Motor Collet

Figure 49: Inserting the Telescoping Tube

Figure 50: Inserting the Dissecting Tool


Figure 51: Locking Rings

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NOTE: MPSS does not provide cranial perforator devices.

Installation:

1. Slide the attachment over the motor collet aligning the triangular markers (Figure 52). You will feel and hear the attachment click into place when it is fully seated.
2. Turn the attachment to the  position (Figure 53).
3. To install a cranial perforator device with a Hudson shank, pull back on the collar of the attachment (Figure 54).
4. Insert the device and release the collar.

NOTE: A cranial perforator device may be installed in the attachment before the attachment is installed on the motor.

Removal:

1. To remove the cranial perforator device, pull back on the collar of the attachment.
2. To remove the attachment, turn it to the  position and slide it off of the motor.

WARNING: Consult the cranial perforator device labeling for the recommended speed specifications.

Maximum Speed of Perforator Attachments

Gas Pressure (Dynamic)	Model AD01 Output Speed (Max.)	Model AD03 Output Speed (Max.)
80 psi	660 rpm	850 rpm
100 psi	820 rpm	1050 rpm
120 psi	890 rpm	1140 rpm

Figure 52: Sliding the Attachment over the Motor Collet


Figure 53: Attachment in the Locked Position


Figure 54: Pulling Back on the Collar

Jacobs® Chuck Attachments

The Jacobs Chuck attachment is a non-cannulated 5/32" chuck with key for drilling.

Installation:

1. Slide the attachment over the motor collet aligning the triangular markers (Figure 55). You will feel and hear the attachment click into place when it is fully seated.
2. Turn the attachment to the  position (Figure 56).
3. To install a drill bit, turn the Jacobs key to open the collar.
4. Insert the drill bit and tighten the collar.

NOTE: A drill bit may be installed in the attachment before the Jacobs Chuck attachment is installed on the motor.

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Figure 55: Sliding the Attachment over the Motor Collet



Figure 56: Attachment in the Locked Position

Bone Mill Attachment

The Legend Bone Mill attachment is composed of a non sterile base and gamma-sterilized, single use disposable components. The gamma-sterilized single use disposable components consist of one bowl, one cap (with LOW and HIGH settings), and one spatula. A Legend or MR7 motor provides the power to drive the Bone Mill attachment. The Bone Mill attachment is intended to mill bone into particles 1 to 5mm in size.

WARNINGS:

- Do not operate the Bone Mill attachment without the bowl, and cap secured in place.
- Do not come in contact with the interior of the disposable bowl and cap during bone milling.
- Do not use the Bone Mill disposable components for more than one surgical procedure, as this may cause cross-contamination and affect patient safety.
- Use only Legend/MR7 motors and Legend Bone Mill disposables in combination with the Bone Mill attachment. Use of other devices will void the manufacturer’s warranty.

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CAUTIONS:

- The addition of any liquid to the milling process should be administered in small volumes to prevent leakage.
 - Do not use Bone Mill disposables to process more than 6 bone pieces approximately 1.5 cc in size at one time.
 - Do not use Bone Mill disposables to process more than a total of 40 cc of bone during any one surgical procedure.
3. Secure the cap onto the bowl by rotating counterclockwise. Do not activate the Bone Mill Attachment without securing the cap.
 4. Use the HIGH cap position when starting with large pieces of bone or large batch amounts. Use the LOW cap position to finish milling a batch or when starting with small pieces or batch amounts (Figure 58).
 5. Insert the motor into the motor connection on the Bone Mill base and rotate the motor to the  position (Figure 59).
 6. Activate the motor to mill bone. Continue to run the Bone Mill Attachment until desired bone particle size is achieved. Adjust motor speed to obtain desired particle size.

Figure 57: Bone Mill Attachment


Figure 58: Bone Mill Cap


Figure 59: Inserting the Motor

NOTE: If bone does not process, deactivate the motor, rotate the motor to the  position, and remove the cap. Reinsert the cap in the HIGH position and activate the motor. If bone spins without milling or jams, divide bone chunks into two batches and process separately.

7. Process the bone for 5–15 seconds and visually check to gauge quality of milling before continuing.
8. Continue to process bone in 5–15 second batches until the majority of particles reach desired size.
9. To remove bone particles stop operation of the motor and rotate it to the  position. Remove the cap by rotating it clockwise. Use the spatula inside the bowl and cap to remove bone particles.

NOTES:

- Milling time is dependent on the density of the patient's bone, the ratio of cancellous and cortical bone, the size of the individual particles and the batch amount of the bone to be processed. Always check bone after 10–15 seconds to gauge quality of milling.
- If the blade becomes inoperable, stop operation of the motor and rotate the motor to the  position. Remove the cap and use the spatula to free bone and/or blades.
- After removing bone particles from the bowl, the milling process may be repeated with additional bone specimens, as needed. (Up to a total of 40 cc of bone may be processed by the Bone Mill Disposable during any one surgical procedure.)

After Use

1. Rotate the motor to the  position.
 2. Unlatch the Bone Mill disposables from the base.
 3. Discard the bowl, cap, and spatula per hospital procedure.
 4. Clean the Bone Mill base according to the directions below.
 5. Clean the motor according to its instruction manual.
-  Bone Mill bowl, cap, and spatula are single use devices. Opened but unused Bone Mill disposables may not be re-sterilized under any sterilization method.

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
Rx Only

Glossary of Universal Symbols


	Attention, see <i>Instructions for Use</i> .
Rx Only	United States federal law restricts this device to sale by or on the order of a physician.
	Reference Number
	Lot Number
	Serial Number
	Quantity
	For single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death
≈	Approximately equal to
	Non-Sterile
	Sterilized by Gamma Irradiation
	Use by date
	Date of manufacture
	Temperature Limitations
	Compliant with European Council Directive MDD 93/42/EEC.
	Authorized Representative in the European Community
	USA Only
	Do not dispose to unsorted municipal waste.
	Manufacturer
	Consult Instructions for Use
	Use with

Table of Contents

TRITON™ Power Surgical Instrument System Components 2

Power Source Requirements 4

“SAFE” Position 5

Pre-Operative Procedures 6

Post-Operative Procedures 8

TRITON™ Power Surgical Instrument System Sterilization 10

TRITON™ High Torque Handpiece 11

TRITON™ High Torque Handpiece Attachments 12

TRITON™ Compact Micro-Saw Handpieces 16

TRITON™ Hose 18

Precautions 19

Warnings 20

Service and Warranty 21

General Warning

Before using the TRITON™ Power Surgical Instrument System, it is EXTREMELY important that you read this instruction manual.

Extreme caution must be observed when using TRITON components near delicate neural and vascular structures. Improper use may cause injury to the patient or hospital staff.

We recommend in-service training by our representatives. If questions arise, please call one of our consultants at 800-643-2773 or 817-788-6400.


Indications for Use

The TRITON™ Power Surgical Instrument System is intended for use during surgical procedures to cut hard tissue, bone, and soft tissue.

TRITON™ Power Surgical Instrument System Components

MIDAS REX® High Speed Motor


- Controlled by foot pedal on TRITON Pneumatic Control.

- 1 – Motor
- 2 – Permanently attached hose
- 3 – Disposable air diffuser
- 4 – Attachment
- 5 – Dissecting Tool


[Refer to Midas Rex Instrumentation System Instruction Manual for detailed operating instructions.]

TRITON™ High Torque (HT) Handpiece

- 1 – “Pistol Grip” Pneumatic Motor
- 2 – TRITON Hose (removable)
- 3 – “Quick Release” button (easy installation and removal of TRITON Attachments)
- 4 – Trigger Control (pressure sensitive to allow variable speed operation)
- 5 – Forward/Reverse Switch (directional)

TRITON™ Compact Micro-Saw Handpieces

- All use TRITON Hose (removable).
- Each has variable speed, Finger Control with a “Safety Slide.”

- 1 – Micro-Oscillating Saw Handpiece
- 2 – Micro-Sagittal Saw Handpiece
- 3 – Micro-Reciprocating Saw Handpiece


TRITON™ Power Surgical Instrument System Components

TRITON™ Pneumatic Control

- Independent operation of both the MIDAS REX® Motor and the TRITON HT Handpiece or a TRITON Compact Micro-Saw Handpiece.
- U.S. version (part #703002) includes internal regulator and operates at 8.3 bar (120 psi) nominal pressure, 13.8 bar (200 psi) maximum pressure.
- International version (part #703003) does not include internal regulator and operates at pressures up to 8.3 bar (120 psi).
- Maximum pressure for the MIDAS REX Motor is either 8.3 bar (120 psi), 10.3 bar (150 psi), or 13.8 bar (200 psi), depending on the pressure requirements of the MIDAS REX Motor being used.
- Also operates MEDNEXT® Drill. Maximum pressure for MEDNEXT Drill is 10.3 bar (150 psi). [Refer to Mednext Bone Dissecting System Instruction Manual (#913001) for detailed operating instructions.]
- Quick-connect couplers for TRITON Auto Lubricator (for use with the MIDAS REX Motor) and for TRITON Hose. Connect MIDAS REX Motor to gold quick-connect coupler on TRITON Auto Lubricator. Use black quick-connect coupler for MEDNEXT Drill.


- 1 – Quick-connect coupler for TRITON Hose
- 2 – Quick-connect coupler for TRITON Auto Lubricator with MIDAS REX Motor or for MEDNEXT Drill

- 3 – DISS (Diameter Index Safety System) fitting for Pneumatic Supply Hose

Pneumatic Supply Hose

DISS fitting connects the gas source to the Pneumatic Control:

- Standard length – 6.1 meters (20 feet)
- Optional length – 9.1 meters (30 feet)

Power Source Requirements

Specifications	TRITON™ HT Handpiece	TRITON™ Compact Micro-Saw Handpiece
Operating pressure	8.3 to 13.8 bar (120 to 200 psi) with U.S. version TRITO Pneumatic Control (#703002); up to 8.3 bar (120 psi) with international version (#703003) 7 bar (100 psi) if connected directly to gas source	8.3 to 13.8 bar (120 to 200 psi) with U.S. version TRITON Pneumatic Control (#703002); up to 8.3 bar (120 psi) with international version (#703003) 7 bar (100 psi) if connected directly to gas source
Activation method	Trigger Control, variable speed	Finger Control, variable speed
Forward/Reverse	Yes	No
Lubrication	None	None
Detachable hose	Yes (blue)	Yes (blue)
Cannulation	3.2 mm (.125")	None

Gas Type

Nitrogen or dry, filtered compressed air.

Required Pressure

U.S. version of TRITON Pneumatic Control with internal regulator (part #703002): 8.3 to 13.8 bar (120 to 200 psi) maximum, 8.3 bar (120 psi) nominal.

International version of TRITON Pneumatic Control without internal regulator (part #703003): 8.3 bar (120 psi) maximum.

Gas Source

Nitrogen cylinder used with a Pressure Regulator (CGA 580 fitting) or in-house compressed air.

"SAFE" Position

TRITON™ High Torque Handpiece

The TRITON HT Handpiece will not operate in the "SAFE" position.


To operate the TRITON HT Handpiece, activate and depress the Trigger Control.

"SAFE": Turn Trigger Control to either side.

"ON": Trigger Control in vertical position will allow activation of handpiece.

"SAFE": Turn directional Forward/Reverse Switch to middle position.

"ON": Forward/Reverse Switch in forward or reverse position will allow activation of handpiece.

TRITON™ Compact Micro-Saw Handpieces

The TRITON Compact Micro-Saw Handpieces will not operate in the "SAFE" position.

To operate Compact Micro-Saw Handpieces, activate and depress the Finger Control.

"SAFE": Lift Finger Control and push "Safety Slide" away from "ON" position.

"ON": To activate, lift Finger Control and push "Safety Slide" to "ON" position.

CAUTION: Always install TRITON Handpieces, Attachments, saw blades and other items with Handpiece in "SAFE" position.

Pre-Operative Procedures


STEP A – Pneumatic Supply Hose Connections (DISS Fittings)

1A. For in-house gas source with DISS fitting, connect one end of Pneumatic Supply Hose to wall fitting.

1B. For in-house gas source with Schrader fitting, connect one end of Pneumatic Supply Hose to DISS-Schrader Adapter. Connect Adapter to wall fitting.

1C. For nitrogen tank gas source, connect one end of the Pneumatic Supply Hose to the Pressure Regulator.


2. Connect the other end of the hose to the Pneumatic Control.

NOTES: When operating TRITON Handpieces without the TRITON Pneumatic Control, connect the male Schrader fitting on TRITON Hose directly to a gas source. Set pressure at 7 bar (100 psi).

When connecting TRITON Pneumatic Supply Hose to a MIDAS REX Safety Valve Regulator (SVR), remove in-line oiler assembly from SVR. Wrench-tighten WF-4 valve to SVR, then hand-tighten Pneumatic Supply Hose DISS fitting to WF-4 valve. Hand-tighten all DISS fittings – DO NOT use wrenches. If DISS fittings are unavailable, contact MIDAS REX® for an adapter.

STEP B – TRITON™ Hose and Handpiece Connections to Pneumatic Control

1. Connect color-coded blue TRITON Hose to blue quick-connect coupler on the Pneumatic Control.


2. Check to see that TRITON HT Handpiece or Compact Micro-Saw Handpiece is in "SAFE" position.


3. Connect TRITON HT Handpiece or Compact Micro-Saw Handpiece to TRITON Hose with "J" lock. Handpieces can be connected while the TRITON Hose is pressurized.

Pre-Operative Procedures


STEP C – MIDAS REX® Connections to Pneumatic Control


1. Unscrew Lubricant Reservoir in "OFF" direction. (If Pneumatic Control is pressurized, turn off gas source and depress foot pedal 3 to 5 times to remove air from system.)

2. Check that siphon tube is seated firmly in side hole under reservoir cap.


3. Fill with MIDAS REX Lubricant until lubricant reaches "FULL" line. Reattach and firmly hand-tighten reservoir. CAUTION: Do not cross-thread reservoir. Use only MIDAS REX Lubricant.

4. Place alignment pin of Auto Lubricator into alignment port of Pneumatic Control. Push Auto Lubricator into black quick-connect coupler.

5. Attach Air Diffuser to base of MIDAS REX motor hose.

6. Connect MIDAS REX Motor to gold quick-connect coupler on the Auto Lubricator.
NOTE:
Connect MEDNEXT® Drill directly to black quick-connect coupler on TRITON™ Pneumatic Control. Refer to the MEDNEXT® Bone Dissecting System Manual for lubrication, cleaning and maintenance of the Mednext Drill System.

7. Pressurize system and adjust lubricant flow rate. Depress Pneumatic Control foot pedal to run motor, then adjust flow control on Auto Lubricator (see chart this page).


LUBRICANT DRIP RATES

- **MIDAS REX Classic Motor**
1 drop per 10-12 seconds
- **MIDAS REX III Motor**
1 drop per 30-60 seconds

Post-Operative Procedures

Disassembly


1. Discard used dissecting tools and blades.

- 2A. Remove MIDAS REX® Motor from Auto Lubricator by pushing down gold quick-connect coupler. Discard Air Diffuser. If Auto Lubricator is stored separately from Pneumatic Control, store upright or remove lubricant from reservoir before storage.

- 2B. Hold TRITON™ Hose firmly and remove by turning blue quick-connect coupler counter-clockwise.


3. Turn off gas source and depress foot pedal 3 to 5 times to remove air from Pneumatic Supply Hose.
4. Remove Pneumatic Supply Hose from Pneumatic Control and gas source.

Post-Operative Procedures


Cleaning

Clean the MIDAS REX® Motor according to instructions in the MIDAS REX Instrumentation System Operative Manual.

CAUTION: A TRITON™ Handpiece must be connected to the TRITON Hose to keep water out of the Handpiece during the cleaning process.

1. Clean TRITON HT Handpiece, TRITON Compact Micro-Saw Handpieces, TRITON Hose and all TRITON Attachments immediately after use with an appropriate mild detergent solution and brush.


2. Rinse instrumentation under running water to remove all traces of detergent solution. If possible, use distilled water for the final rinse.

3. Disconnect the TRITON Hose from the TRITON HT Handpiece or Compact Micro-Saw Handpiece. Dry all components with a lint-free towel.

- 1A. Clean central cannulation shaft of the TRITON HT Handpiece with a wire driver cannulation brush. Clean nose of the TRITON HT Handpiece, Reciprocating Saw Attachment and Compact Micro-Reciprocating Saw Handpiece with the same brush.

- 1B. Flush the nose of the TRITON HT Handpiece, Reciprocating Saw Attachment and Compact Micro-Reciprocating Saw Handpiece with a water-pik.

WARNINGS: DO NOT lubricate or oil the TRITON Handpieces or the TRITON Hose. DO NOT immerse. DO NOT use ultrasonic washer/sterilizer.


If any TRITON Handpiece is ACCIDENTALLY IMMERSED in saline, disinfectant, cleaning fluid, or any other corrosive substance, take the following steps:

1. Immerse the Handpiece in distilled water for 1 minute to dilute the corrosive fluid.
2. Immediately, steam sterilize the Handpiece in a prevacuum sterilizer at 132° C for 4 minutes.


TRITON™ Power Surgical Instrument System Sterilization

Sterilization Guidelines

1. Remove TRITON™ Hose from TRITON™ Handpiece.
DO NOT attempt to remove hose from the MIDAS REX® Motor or MEDNEXT® Drill.

2. Place cleaned MIDAS REX® Motor or MEDNEXT® Drill and TRITON™ Hose in bottom of instrument case.

3. Place TRITON™ Handpieces, attachments, and accessories in appropriate depressions of insert trays.
Set tray in case.

NOTE:

Midas Rex or Mednext attachments can be placed in their respective insert trays and positioned in the TRITON instrument case.

4. Place lid on instrument case and sterilize.

Exposure Times

Pre-Vacuum Steam Sterilization
132°C/270°F for 4 minutes.

Gravity Displacement Steam Sterilization
132°C/270°F for 30 minutes.

WARNINGS:


Following steam sterilization, the instrumentation must be allowed to cool to room temperature prior to use.

Dissecting tools and saw blades are single use, disposable items.


DO NOT steam sterilize Pneumatic Control, Auto Lubricator, Pneumatic Supply Hose, DISS to Schrader Adapter or Pressure Regulator; the sterilization process could damage components and void the warranty.

TRITON™ High Torque Handpiece


TRITON™ High Torque Handpiece Operation


1. Attachments may be preloaded before insertion into the TRITON HT Handpiece.

2. With TRITON HT Handpiece in "SAFE" position, insert preloaded Attachment by depressing quick-release button on top of Handpiece. Snap Attachment into Handpiece with a slight twisting motion until seated.


3. Place the TRITON HT Handpiece in "RUN" position with the Trigger Control vertical. The Forward/Reverse Switch should be in the "FORWARD" position.

4. After use, return Trigger Control to "SAFE" position prior to removing Attachment. Turn Forward/Reverse Switch to middle position.

5. Remove Attachment by depressing quick-release button on top of TRITON HT Handpiece.


NOTE: The TRITON HT Handpiece has a Handle Extension that screws into the back of the Handpiece. The Handle Extension provides balance and two-handed control for various drilling and cutting applications.


CAUTION: Insert all Attachments into TRITON HT Handpiece with Handpiece in "SAFE" position.

TRITON™ High Torque Handpiece Attachments


Sagittal Saw Attachment


1. Insert Sagittal Saw Attachment into TRITON HT Handpiece in a position to allow easy insertion of the Sagittal Saw Key. Note that Attachment may be installed in 12 different positions to facilitate proper surgical access.

2. Insert the Sagittal Saw Key into the Attachment and turn counterclockwise until there is slight resistance.

3. Insert the blade in the space between the two jaws, ensuring that the blade is fully seated.

4. Turn the Sagittal Saw Key clockwise to lock the blade. Run briefly, then retighten blade.

CAUTION: Do not over-tighten.


To remove the blade, insert the Sagittal Saw Key into the Attachment and turn counterclockwise.

TRITON™ High Torque Handpiece Attachments


Reciprocating Saw Attachment

1. Insert Reciprocating Saw Attachment into TRITON HT Handpiece while Handpiece is in "SAFE" position. The Attachment may be installed in different positions to facilitate proper surgical access.

2. With collet nut loosened, insert blade until fully seated.


3. Finger-tighten collet nut. Run briefly, then retighten collet nut.


To remove the blade, unscrew the collet nut.

TRITON™ High Torque Handpiece Attachments

AO/Synthes® Chuck and Trinkle Chuck Attachments


AO/Synthes® Chuck (left) and Trinkle Chuck (right) Attachments.

To install drill bit, pull back on Attachment Collar. Insert drill bit and release Attachment Collar. Reverse procedure to remove drill bit.


Jacobs Chuck Attachments

Jacobs Chucks and key placement for 5/32" and 1/4" HT sizes.


To install drill bit, turn key to open ridged collar. Insert drill bit. Reverse procedure to remove drill bit.


Keyless Jacobs Chucks are available in 1/8" and 1/4" HT sizes.

Hudson and Zimmer® Chuck HT Attachments

Hudson Chuck HT (left) and Zimmer® Chuck HT (right) Attachments.

To install instrument, pull back on Attachment Collar, then insert male end of instrument into Chuck. Reverse procedure to remove drill bit.

TRITON™ High Torque Handpiece Attachments

Wire and Pin Collet Attachments

Wire Collet accepts wires up to 1.6mm (.062") in diameter. Pin Collet accepts pins up to 3.2mm (.125") in diameter.


1. Insert Wire or Pin Collet while TRITON HT Handpiece is in "SAFE" position.


2. Screw the Cannulated Extension into back of the TRITON HT Handpiece to protect operator from the point of the wire or pin.

3. Insert wire or pin into the front or back of the TRITON HT Handpiece.

4. Put the instrument in the "RUN" position by positioning the Trigger Control vertically. Turn Forward/Reverse Switch to the "FORWARD" position.

5. Squeeze wire/pin Advance Lever and hold it down.

6. Depress Trigger Control to drive the wire/pin. The pressure-sensitive trigger allows variable speed operation.
7. To obtain additional wire/pin length, release the wire/pin Advance Lever and Trigger Control. Pull back on the instrument. Squeeze the wire/pin Advance Lever and Trigger Control to drive the wire.
8. To remove threaded wire/pin, put the Forward/Reverse Switch in "REVERSE," squeeze wire/pin Advance Lever and depress Trigger Control.
9. Remove Wire or Pin Collet while TRITON HT Handpiece is in "SAFE" position.


TRITON™ Compact Micro-Saw Handpieces

Compact Micro-Sagittal Saw Handpiece

Insert blade into Compact Micro-Sagittal Saw Handpiece with Safety Slide in the “SAFE” position.

1. To install blade, lift Saw Blade Lever.

2. Align hole in blade with alignment pin in vice.


3. Push Saw Blade Lever down to lock blade in place. DO NOT force Saw Blade Lever.

4. To activate, depress Finger Control with Safety Slide in “ON” position.
5. To remove blade, push Safety Slide into “SAFE” position and lift Saw Blade Lever.


CAUTION: Use Micro-Sagittal Saw Blades only. DO NOT USE HT-SAGITTAL SAW BLADES.

Compact Micro-Reciprocating Saw Handpiece*

Insert blade into Compact Micro-Reciprocating Saw Handpiece with Safety Slide in the “SAFE” position.

1. With collet nut loosened, insert blade until fully seated.

2. Finger-tighten collet nut. Run briefly, then retighten collet nut.

3. To activate, depress Finger Control with Safety Slide in “ON” position.
4. To remove blade, push Safety Slide into “SAFE” position and unscrew collet nut.

Compact Micro-Oscillating Saw Handpiece

Insert blade into Compact Micro-Oscillating Saw Handpiece with Safety Slide in the “SAFE” position.


1. To install blade, use Oscillating Key to unscrew hex nut.

2. Insert blade underneath hex nut and washer. Tighten hex nut with Oscillating Key. Run saw briefly, then retighten hex nut.

3. To activate Micro-Oscillating Saw Handpiece, depress Finger Control with Safety Slide in “ON” position.

4. To remove blade, push Safety Slide into “SAFE” position and use Oscillating Key to unscrew hex nut.

CAUTION: When connecting any Compact Micro-Saw Handpiece to the TRITON Hose, ensure the Safety Slide is in the “SAFE” position. Insert blades into, or remove blades from, any Compact Micro-Saw Handpiece with Safety Slide in the “SAFE” position.

TRITON™ Hose

TRITON™ Hose:

1. **Inspect TRITON Hose before each use:** Check carefully for any punctures, tears, or worn areas. Also check the metal bands at each end of the hose for looseness or wear. Do not use a hose that has any of these conditions.
WARNING: Do not attempt to repair a torn hose by wrapping tape around it. This could restrict the exhaust flow and cause the hose to explode.
2. **Test the TRITON Hose prior to surgery:**
 - Wear safety goggles or glasses.
 - Attach the TRITON Hose to the Pneumatic Control.
 - Slowly open the nitrogen control valve (nitrogen tank with regulator or “in-house” nitrogen). Listen for any sounds of leakage. Leakage indicates a worn or damaged hose. **Do not use hose if there is any evidence of pressure leakage.**
 - Attach the TRITON Handpiece to the hose and test for normal operation. Note: if the air hose does not easily insert into the handpiece there are two possible problems: 1) The hose fitting or the handpiece fitting may be bent or damaged. Try another hose and handpiece to see if the damaged part can be identified. It must be returned for service. 2) The O-ring in the handpiece connector may be dried out. It can be lubricated with a small amount of STERILE petroleum jelly or sterile water (not saline).
 - Ensure that the hose does not slip loose from the handpiece when in use. It could whip about, possibly causing an injury.
3. **Return the TRITON Hose at least once per year for routine inspection and service.** Aging and repeated usage/sterilization will gradually deteriorate the hose materials. Eventually, the materials may crack and break under pressure. This could result in the hose exploding and whipping about, possibly causing an injury. This is why it is important to routinely return the hose for inspection and service.

Precautions

- Prior to each use, the entire TRITON™ Power Surgical Instrument System must be inspected for proper operation. Check all components for gas leakage and return for servicing if necessary.
- DO NOT OPERATE power system if the TRITON Hose is damaged.
- Inspect the TRITON Hose before each use for any punctures, tears, or worn areas. Also check the metal bands at each end of the hose for looseness or wear. Do not attempt to repair a torn hose by wrapping tape around it; such taping could restrict the exhaust flow and cause the hose to rupture. If the TRITON Hose is damaged, send it in for repair.
- Return the TRITON Hose at least once per year for routine inspection and service. Normal aging, repeated usage, and sterilization will gradually deteriorate the hose materials, causing cracks and breakage under pressure.
- Attach the TRITON Handpieces to the hose and test for normal operation. NOTE: If the TRITON Hose does not insert easily into the Handpieces, there are two possible problems:
 1. The hose fitting or the Handpiece fitting may be bent or damaged. Try another hose and Handpiece to see if the damaged part can be identified. Return damaged parts for service.
 2. The o-ring in the Handpiece connector may be dried out. It can be lubricated with a small amount of sterile petroleum jelly or sterile water (not saline).
- Wear safety goggles or glasses when using the TRITON Power Surgical Instrument System.
- DO NOT initiate patient contact without testing instruments completely. Test for vibration, heat, etc. in the Handpieces, Attachments, and hoses.
- ALWAYS change Handpieces and Attachments with instrument in the "SAFE" position.
- DO NOT operate Saw Attachments or Compact Micro-Saws if blades cannot be firmly secured.
- DO NOT use excessive force during surgery. Using excessive force may cause vibration or breakage, which may result in serious injury.
- USE ONLY TRITON saw blades with the TRITON HT Handpiece Saw Attachments and Compact Micro-Saw Handpieces. Exact fit of blades is essential for stability and proper operation of equipment.

Warnings

Warnings

- **DO NOT** operate the TRITON™ Power Surgical Instrument System unless you have been trained by an official Medtronic Midas Rex representative.
- Prior to each use, the entire TRITON Power Surgical Instrument System must be inspected for proper operation. Check all components for air leakage and return for servicing if leakage is noticed.
- **DO NOT** initiate patient contact without testing instruments completely: test for vibration, heat, etc. in the handpiece, hose, and attachments.
- Eye protection should be worn by operating room personnel and surgeons during use of the TRITON Power Surgical Instrument System.
- Following steam sterilization, the components must be allowed to cool to room temperature prior to use.
- **DO NOT** operate power system if TRITON Hose is damaged.
- **DO NOT** operate handpiece without making sure the proper attachment is used.
- **DO NOT** use an attachment that gets hot or an attachment that is damaged in any way.
- **DO NOT** operate the Sagittal Saw Attachment, Reciprocating Saw Attachment, Compact Micro-Sagittal Saw, Compact Micro-Oscillating Saw, or Compact Micro-Reciprocating Saw if blades cannot be firmly secured.
- **DO NOT** use excessive force during surgery. Excessive force may lead to vibration or breakage which may result in serious injury.
- Use **ONLY** blades manufactured by Medtronic Midas Rex. Exact fit of blades is essential for stability and proper operation of equipment.
- Always install a new blade prior to use.

!USA **Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

Contact: Medtronic Midas Rex
Technical Support: 800-643-2773
817-788-6400

CAUTION ABOUT REPAIRS

When MIDAS REX®, MEDNEXT® or TRITON™ equipment requires servicing or repair, it is strongly urged that the equipment be returned to MEDTRONIC MIDAS REX, Fort Worth, TX, for quality assured service by factory-trained personnel who will utilize genuine parts as required.

LIMITED WARRANTY

Medtronic Midas Rex warrants that the instruments and parts sold by Medtronic Midas Rex have been tested, inspected and shipped from the factory in proper working condition, free of visible defects.

Medtronic Midas Rex warrants the following to be free from defects in material and workmanship for a period of one (1) year after the date of purchase:

MIDAS REX® Motors*	TRITON™ Pneumatic Controls
MIDAS REX® Hoses	TRITON™ Attachments
MIDAS REX® Foot Controls	Bur Storage Cabinets
MEDNEXT® Handpieces	Bur Racks
MEDNEXT® Hoses	Pressure Regulators
MEDNEXT® Foot Controls	Sterilization Trays
TRITON™ Handpieces	Pneumatic Supply Hoses
TRITON™ Hoses	

Any service or repair work done by Medtronic Midas Rex on the equipment listed above, following the one (1) year warranty period, will be free from defects in materials and workmanship for a period of six (6) months after the date of shipment back to the customer.

*MIDAS REX® Driver Motor

Medtronic Midas Rex warrants the MIDAS REX® Driver Motor (Product Number: MRDM) to be free from defects in material and workmanship for a period of six (6) months after the date of purchase.

Any service or repair work done by Medtronic Midas Rex on the MIDAS REX® Driver Motor (MRDM), following the six (6) month warranty period, will be free from defects in materials and workmanship for a period of ninety (90) days after the date of shipment back to the customer.

Medtronic Midas Rex warrants the following to be free from defects in material and workmanship for a period of ninety (90) days after the date of purchase:

- MIDAS REX® Attachments and Tubes
- MEDNEXT® Attachments

Medtronic Midas Rex warrants all other accessories not listed above to be free from any defects in material or workmanship for a period of ninety (90) days after the date of purchase. Medtronic Midas Rex will service or replace any accessories found to be defective on receiving the defective part.

Service and Warranty (continued)

This warranty does not apply to any instrument or part, which has been:

1. Damaged by accident, abuse, misuse, or tampering.
2. Modified without written permission of Medtronic Midas Rex.
3. Used contrary to the instructions and warnings specified in the Medtronic Midas Rex instruction manuals for MIDAS REX®, MEDNEXT® or TRITON™ equipment.
4. Serviced or repaired by anyone other than Medtronic Midas Rex.
5. Damaged by the use of accessories other than genuine MIDAS REX® Tools and Attachments, MEDNEXT® Burs and Attachments, or TRITON™ Blades and Attachments.

This warranty does not include any implied warranties or any warranty of any type including without limitation merchantability or fitness for a particular purpose, other than those expressly contained herein. This warranty is exclusive and in lieu of other warranties, whether written, oral or implied.

LIABILITY IN CLINICAL APPLICATIONS

While Medtronic Midas Rex guarantees complete compatibility among its products within a specific product line, it disclaims any responsibility should a surgical procedure be complicated by a misadventure if:

1. The MIDAS REX® Motors are used with dissecting tools other than MIDAS REX® Dissecting Tools that are manufactured, authorized, or quality controlled by the Medtronic Midas Rex.
2. MIDAS REX® Dissecting Tools manufactured, authorized, or quality controlled by the Medtronic Midas Rex are used with any Motor other than the appropriate MIDAS REX® motor.
3. The MEDNEXT® Drills are used with burs other than MEDNEXT® Burs that are manufactured, authorized, or quality controlled by the Medtronic Midas Rex.
4. MEDNEXT® Burs manufactured, authorized, or quality controlled the Medtronic Midas Rex are used with any drill other than the appropriate MEDNEXT® Drill.
5. The TRITON™ HT Handpiece is used with attachments or saw blades other than TRITON™ attachments or saw blades that are manufactured, authorized, or quality controlled by the Medtronic Midas Rex.
6. The TRITON™ Compact Micro-Saws are used with saw blades other than TRITON™ saw blades that are manufactured, authorized, or quality controlled by the Medtronic Midas Rex.
7. TRITON™ saw blades manufactured, authorized, or quality controlled by the Medtronic Midas Rex are used with any handpiece/attachment other than the appropriate TRITON™ handpiece/attachment.

Note: The Limited Warranty does not apply to services required as the result of damage caused by the use of non-genuine MIDAS REX®, MEDNEXT® or TRITON™ manufactured cutting tools. Peak performance, reliability and maximum service life from your MIDAS REX®, MEDNEXT® or TRITON™ instrumentation may be assured by using only those products that are manufactured by and sold through Medtronic Midas Rex, Fort Worth, Texas.

MANUFACTURER'S RECOMMENDATION

It is strongly recommended that Surgeons and Operating Room Personnel become thoroughly familiar with all Medtronic Midas Rex equipment prior to use in surgery. To assist in this effort, Medtronic Midas Rex offers several Hands-On Workshops. Contact your local Medtronic Neurologic Technologies Sales Representative or Medtronic Midas Rex for schedules.

©2011 Medtronic, Inc.
All rights reserved
Printed in the USA

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Medtronic Powered Surgical Solutions
4620 North Beach Street
Fort Worth, Texas 76137
Telephone: 800-433-7080
817-788-6488

Authorized Representative
in Europe:
Medtronic B. V.
Earl Bakkenstraat 10
6422 PJ Heerlen
THE NETHERLANDS

- ▣ See instructions for use.
- ▣ Vedere istruzioni per l'uso
- ▣ Ver instrucciones para el uso
- ▣ Bedienungsanleitung beachten
- ▣ Voir notice d'utilisation

ENGLISH**IMPORTANT INFORMATION FOR TRITON™ SURGICAL INSTRUMENTS****PURPOSE:**

This instrument is intended for use in surgical procedures.

WARNING NOTICE:

Review the instruction manual for the TRITON power system before operating the instrument. Failure to do so could result in injury to the patient or operating room staff.

If there is any doubt or uncertainty concerning the proper use of this device, please contact MEDTRONIC POWERED SURGICAL SOLUTIONS Customer Service for instructions.

GENERAL WARNINGS:

1. This device should not be used for purposes other than those for which it is recommended.
2. This device must not be used for alternative purposes, since this may overstress or otherwise damage it such that it may not be suitable to perform its designed function. Do not use this device for any action for which it was not intended, such as hammering, prying, or lifting.
3. To avoid injury, this device should be carefully examined prior to use for functionality or damage. A damaged device should not be used. Additional backup devices should be available in case of unexpected need.
4. If this device is improperly loaded into the driving instrument, it may slip distally and cause patient injury.
5. Eye protection should be worn by operating personnel using this device.

POSSIBLE ADVERSE EFFECTS:

Breakage, slippage, misuse, or mishandling of this device may cause injury to the patient or operating personnel.

Improper maintenance or handling can render the device unsuitable for the intended purpose, or even dangerous to the patient or operating personnel.

Extreme care should be taken to ensure that this device remains in good working order. Misuse of the device causing corrosion, bending, and/or fracture of any sections of the device may inhibit or prevent proper function.

It is important that the surgeon exercise extreme caution when working in close proximity to vital organs, nerves, or vessels, and that the force applied when correcting the position of the device is not excessive, such that it might cause injury to the patient.

CAUTION: FOR USE BY OR ON THE ORDER OF A PHYSICIAN.

CAUTION: FEDERAL (US) LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN ONLY.

This device should be used only by physicians familiar with the device, its intended use, any additional instrumentation and any available surgical techniques.

This device should be used only in conjunction with other TRITON products.

PACKAGING, SHIPPING, AND STORAGE:

The packaging should be intact upon receipt. Devices in damaged packaging should not be used. Do not ship or store with other devices that may crush or damage the product. Devices should be shipped in containers such as corrugated boxes that prevent damage. Devices should be stored in a dry environment, protected from direct sunlight.

EXAMINATION:

This device must be examined by the user prior to use in surgery.

Examination should be thorough, and in particular, should take into account visual and functional inspection of the working surfaces and the presence of any cracks, bending, or distortion and that all components of the device are complete.

Never use this device with obvious signs of excessive wear, damage or that are incomplete or otherwise non-functional.

FURTHER INFORMATION:

No warranties expressed or implied are made. No claims of merchantability, fitness, or use other than those expressly contained herein and in the TRITON instruction manual. If further information on the use of this device is needed, please check with MEDTRONIC POWERED SURGICAL SOLUTIONS Customer Service at the addresses below.

PRODUCT COMPLAINTS:

Any Health Care Professionals (e.g. customer users of MEDTRONIC POWERED SURGICAL SOLUTIONS instruments) who have any complaint or who have experienced dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, and/or performance should notify the distributor, MEDTRONIC POWERED SURGICAL SOLUTIONS. Furthermore, if any instrument does not meet any of its performance specifications or otherwise does not perform as intended or is suspected of doing so, the distributor, MEDTRONIC POWERED SURGICAL SOLUTIONS, should be notified immediately. If any MEDTRONIC POWERED SURGICAL SOLUTIONS product ever may have caused or contributed to the death or serious injury of a patient, the distributor, MEDTRONIC POWERED SURGICAL SOLUTIONS, should be notified as soon as possible by telephone or written correspondence. When filing a complaint, please provide the component(s) name and number, lot or serial number, your name and address, and the nature of the complaint.

COMPLAINTS SHOULD BE ADDRESSED TO:

Regulatory Affairs
Medtronic Powered Surgical Solutions Telephone: 800-643-2773 or
4620 North Beach Street 817-788-6400
Fort Worth, Texas 76137
USA

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
9 cm attachment length for use in 9-M attachment
Match Head dissecting tool tip design

3.0 mm dissecting tool head diameter
Classic legacy nomenclature (reference only)
7 cm attachment length for use in 7-6ST straight or angled attachment
Ball dissecting tool tip design
Mednext tip design

4.0 mm dissecting tool head diameter
Mednext legacy nomenclature (reference only)
8 cm attachment length for use in F2-B1 footed attachment or 8-B attachment
2.3 mm dissecting tool head diameter
Classic legacy nomenclature (reference only)
8 cm attachment length for use in 8-B attachment
Twist Drill dissecting tool tip design
4.0 mm depth stop
1.1 mm dissecting tool head diameter

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Medtronic

Medtronic

Medtronic Powered Surgical Solutions (MPSS)

4620 N. Beach Street

Fort Worth, Texas 76137


USA

www.medtronic.ru

tel 817.788.6400

fax 817.788.6225

Custodian:

Jenna Groves

Senior Regulatory Affairs Specialist

Medtronic Powered Surgical Solutions

28-July-2017

Date


Additional Instruction Manual**Surgical Pneumatic Instruments with Accessories**ATTESTATION

STATE OF TEXAS

TARRANT COUNTY


On this 28 day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, **Additional Instruction Manual**, issued by Medtronic Powered Surgical Solutions, presented to me by the document's custodian, **Jenna Groves**, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
- TRITON Sterilization Tray, model 801102 – not more than 20 pcs.
 - Midas Rex Legend 1.2mm Cleaning Brush, models: PA300, PA310, PA320 – not more than 100 pcs.
5. Midas Rex MR7 Touch:
- Midas Rex MR7 Touch Motor, model: PM710
 - Midas Rex Legend Lubricant/Diffuser Pack, model: PA700 – not more than 300 pcs.
 - Pneumatic Control Unit, models: PC100, PC700 – not more than 20 pcs.
 - Legend Regulator DISS conn, model PA200 – not more than 20 pcs.
 - Midas Rex Legend MR7 Regulator Hose, model: PC710 – not more than 20 pcs.
 - Midas Rex Motor Hoses, models: PC230, PC231 – not more than 20 pcs.
 - Midas Rex Adapters, models: PC110, PC120, PC140, PC720 – not more than 20 pcs.
 - Midas Rex WF4 Connector, model: PC130.
 - Operational documentation.
- Midas Rex MR7 Touch Accessories
- Midas Rex Motor Wrench, model: PA130 – not more than 50 pcs.
 - Midas Rex Legend 3/4 DIN Instrument Kit, models: PA400, PA600, PA500 – not more than 20 pcs.
 - TRITON Sterilization Tray, model 801102 – not more than 20 pcs.
 - Midas Rex Legend 1.2mm Cleaning Brush, models: PA300, PA310, PA320 – not more than 100 pcs.
6. TRITON:
- TRITON Motor high Torque, model: 720001
 - Pneumatic Control Unit, models: PC100, PC700 – not more than 20 pcs.
 - Legend Regulator DISS conn, model PA200 – not more than 20 pcs.
 - Midas Rex Legend MR7 Regulator Hose, model: PC710 – not more than 20 pcs.
 - Triton hand piece hose, models 707020, 720101 – not more than 20 pcs.
 - Midas Rex Adapters, models: PC110, PC120, PC140, PC720 – not more than 20 pcs.
 - Midas Rex Motor Hoses, models: PC230, PC231 – not more than 20 pcs.
 - Midas Rex WF4 Connector, model: PC130.
 - Operational documentation.
- TRITON Accessories:
- Midas Rex Motor Wrench, model: PA130 – not more than 50 pcs.
 - Midas Rex Legend 3/4 DIN Instrument Kit, models: PA400, PA600, PA500 – not more than 20 pcs.

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


схема/Design and Technical
scheme

INDEX	QTY.	PART NUMBER	DESCRIPTION
1	1	100-075-20	MOTOR SEAL 2
2	1	100-510-05	ROTOR SHAFT-3 VANE OFFSET, LEGEND GOLD, RULON J SEAL
3	1	100-510-08	ROTOR HOUSING, LEGEND GOLD MOTOR
4	1	100-510-17	FERRULE, EXHAUST, MOTOR A
5	1	100-510-18	DISTRIBUTOR, MOTOR A
6	3	100-510-31	VANE, LEGEND GOLD MOTOR
7	1	100-510-38	REAR BEARING HOUSING, LEGEND GOLD MOTOR
8	1	100-510-50	SWIVEL CONVERTIBLE, HOSE END, MOTOR A
9	1	100-510-54	FRONT SWIVEL UNION, LEGEND GOLD MOTOR
10	1	100-510-55	COMPRESSION SHIM, MOTOR A
11	1	100-510-56	NUT, COMPRESSION, LEGEND GOLD MOTOR
12	1	100-510-57	NOZZLE, FRONT, MOTOR A
13	1	100-510-63	SOUND TUBE SPACER, GOLD MOTOR
14	1	100-510-64	HOUSING, FRONT BEARING, LEGEND GOLD MOTOR
15	1	100-510-66	BARB, HP HOSE, GOLD, CRIMPED
16	1	100-550-1	CASE, MOTOR, LEGEND GOLD
17	1	N710423B001	ASSEMBLY, DOUBLE O COLLET, PNEUMATIC
18	1	120-030-1	OIL, MPSS REX LUBRICANT
19	1	120-051-3	EXHAUST TUBE, OILER
20	1	120-051-17	I.D LONG EXHAUST HOSE, LEGEND PLATINUM
21	1	N708764B001	ASSEMBLY, HOSE, HIGH PRESSURE
22	1	120-051-23	SPRING, COMPRESSION, MODIFIED, .430 OD X .041 WIRE X 10" LONG
23	1	120-051-24	SEAL, LIP, OILER HOUSING
24	1	120-051-25	STEM, LUBRICATOR/DIFFUSER, LEGEND
25	1	120-051-26	SCREEN, LUBRICATOR/DIFFUSER, LEGEND
26	1	120-051-28	CONNECTOR SPACER, OILER HOUSING, LEGEND
27	1	N405284B001	EXHAUST HOSE, LEGEND
28	1	120-051-30	PRESSURE TUBE, OILER, CRIMPED
29	1	120-051-33	HOUSING, OILER
30	1	120-057	ASSEMBLY, CAPILLARY TUBE, LEGEND GOLD

31	1	900-188	BEARING, .250 ID X .500 OD
32	1	900-2-5	BEARING, .125 ID X .3125 OD
33	1	901-001	BALL, .078 DIA, 440C, SS
34	56	901-009	BALL, 3/64" (.0469) DIAMETER
35	1	902-013	PIN, SPIRAL, .049 X .500
36	1	903-012	O-RING, .063 CS, 0.799 ID (FLUOROELASTOMER)
37	1	903-016	O-RING, .031 K-SEC X .312 ID, 70 DURO SILICONE
38	1	903-017	O-RING, .032 K-SEC X .062 ID, 70 DURO SILICONE
39	1	903-030	O-RING, .079 CS, X .945 I.D.
40	1	903-037	O-RING, .032 K-SEC, .220 I.D., 70 DURO SILICONE
41	1	903-038	O-RING, .040 K-SEC, .226 I.D., 70 DURO SILICONE
42	1	903-039	O-RING, .040 K-SEC, .226 I.D., 70 DURO SILICONE
43	1	903-040	O-RING, .051 K-SEC, .547 I.D., 70 DURO SILICONE
44	1	903-30-622	O-RING, .030 CS X .622 ID
45	1	903-086	SPRING, .062 O.D. X .015 WIRE X .160 LONG
46	1	906-050	PLUG, 1/4" MALE BSP
47	2	906-063	FERRULE, BRASS, .312 ID X .50 L X .225 HOLE, HPH
48	1	910-001	LOCTITE 600
49	3	910-004	WIRE, STAINLESS STEEL, .020 DIAMETER
50	1	910-014	LOCTITE 609
51	1	910-015	PRIMER "N"
52	1	910-017	SUPER LUBE SYNTHETIC GREASE W/PTFE
53	1	910-096	DC 732 RTV SILICONE ADHESIVE SEALANT (NOT SHOWN)
54	1	910-098	LOCTITE 545 THREAD SEALANT
55	1	910-120	SLEEVE, NOISE REDUCTION, LEGEND MOTORS
56	1	911-003	SET SCREW, 2-56 X .003
57	2	911-006	SCREW, SET, SOCKET HEAD, 2-56 X 1/16"
58	2	911-041	SET SCREW, SOCKET HEAD, 4-40 X 1/16"
59	3	912-036	SCREW, 4-40 X 3/16 SS SOCKET HEAD CAP
60	1	913-004	SCREW, PAN HEAD, 2-56 X .125

1. TITLE 2. DATE 3. BY 4. CHECKED 5. APPROVED 6. REVISIONS 7. COMMENTS	8. REV 9. DATE 10. BY 11. CHECKED 12. APPROVED 13. REVISIONS 14. COMMENTS	15. TITLE 16. DATE 17. BY 18. CHECKED 19. APPROVED 20. REVISIONS 21. COMMENTS	22. TITLE 23. DATE 24. BY 25. CHECKED 26. APPROVED 27. REVISIONS 28. COMMENTS
--	---	---	---

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
NOTES:

- 1) ENSURE PARTS ARE FREE OF BURRS AND FIT/FUNCTION WELL. CLEAN ALL METAL COMPONENTS TO REMOVE ANY OIL, DIRT, GREASE OR DEBRIS BEFORE ASSEMBLY AS NEEDED.
- 2) ASSEMBLE PER WORK INSTRUCTIONS SPECIFIED IN ACCOMPANYING ROUTER.
- 3) APPLY MIDAS REX LUBRICANT (ITEM 18) TO O-RINGS PRIOR TO INSTALLATION.
- 4) APPLY PRIMER (ITEM 51) TO THREAD AND ALLOW TO DRY. APPLY ADHESIVE (ITEM 50) AND INSTALL.
- 5) APPLY THREAD SEALANT (ITEM 54) TO MALE THREADS OF SPACER (ITEM 26) AND INSTALL TO OILER HOUSING (ITEM 29).
- 6) PRESS SPACER (ITEM 13) ONTO HOSE END (ITEM 8).
- 7) BRAID SHALL BE WITHIN 0.15" FROM HOSE END BEFORE CRIMPING
- 8) FERRULE SHALL BE PRESSED FLUSH AGAINST FLAT ON BARB OR PRESSURE TUBE.
- 9) COAT WITH SILICONE ADHESIVE (ITEM 53) TO PROTECT FROM WIRE ENDS AND SECURE HOSE.
- 10) LUBRICATE SEAL (ITEM 1) WITH LUBE (ITEM 52) PRIOR TO INSTALLING ROTORSHAFT (ITEM 2).
- 11) REFERENCE 600-023 FOR DEFINITION AND DELINIATION OF MIDAS TRUE TYPE FONT.
- 12) APPLY PRIMER (ITEM 51) TO THREAD AND ALLOW TO DRY. APPLY ADHESIVE (ITEM 48) AND INSTALL.
- 13) QUANTITY BASED ON TYPICAL DIMENSIONS. ADJUST FOR BEST PERFORMANCE.
- 14) LOCATE APPROXIMATELY AS SHOWN.


Midtronics			
SEE RELEVANT PARTS SPECIFICATIONS			
REV. 1.0	REV. 2.0	REV. 3.0	REV. 4.0
TITLE			
ASSEMBLY, FINAL LEGEND GOLD MOTOR			
DRAWING NUMBER			
2175-0001			
REV			
P			

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
SEE DETAIL D

SECURE HOSE (ITEM 27) TO SLEEVE (ITEM 55)

INSTALL OVER SAFETY WIRE SECURING HOSE TO OILER

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

NOTES:


- 1) ENSURE PARTS ARE FREE OF BURRS AND FIT/FUNCTION WELL. CLEAN ALL METAL COMPONENTS TO REMOVE ANY OIL, DIRT, GREASE OR DEBRIS BEFORE ASSEMBLY AS NEEDED.
- 2) ASSEMBLE PER WORK INSTRUCTIONS SPECIFIED IN ACCOMPANYING ROUTER.
- 3) APPLY MIDAS REX LUBRICANT (ITEM 22) TO O-RINGS PRIOR TO INSTALLATION.
- 4) APPLY PRIMER (ITEM 60) TO THREAD AND ALLOW TO DRY. APPLY ADHESIVE (ITEM 59) AND INSTALL.
- 5) APPLY THREAD SEALANT (ITEM 63) TO MALE THREADS OF SPACER (ITEM 30) AND INSTALL TO OILER HOUSING (ITEM 33).
- 6) PRESS SPACER (ITEM 13) ONTO HOSE END (ITEM 11).
- 7) BRAID SHALL BE WITHIN 0.15" FROM HOSE END BEFORE CRIMPING
- 8) FERRULE SHALL BE PRESSED FLUSH AGAINST FLAT ON BARB OR PRESSURE TUBE.
- 9) COAT WITH SILICONE ADHESIVE (ITEM 62) TO PROTECT FROM WIRE ENDS AND SECURE HOSE.
- 10) LUBRICATE SEAL (ITEM 1) WITH LUBE (ITEM 61) PRIOR TO INSTALLING ROTORSHAFT (ITEM 2).
- 11) REFERENCE 600-023 FOR DEFINITION AND DELINIATION OF MIDAS TRUE TYPE FONT.
- 12) APPLY ADHESIVE (ITEM 62) BETWEEN FINGER CONTROL LEVER MOUNT (ITEM 8) AND MOTOR CASE (ITEM 20) TO FORM AN AIRTIGHT SEAL. CLEAN RESIDUE FROM CASE.
- 13) APPLY PRIMER (ITEM 60) TO THREAD AND ALLOW TO DRY. APPLY ADHESIVE (ITEM 57) AND INSTALL.
- 14) QUANTITY BASED ON TYPICAL DIMENSIONS. ADJUST FOR BEST PERFORMANCE.


MOTOR MARKING DETAIL
SCALE 3:1

MARK WITH 'CE MARK', KEYSTROKE "e", MIDAS TRUE TYPE FONT, .3" SQUARE, CENTERED AS SHOWN.

CE MARK DETAIL
SCALE 2:1

PRIMER (ITEM 60) THREADS OF COMPRESSION NUT. APPLY ADHESIVE (ITEM 59). ASSEMBLE COLLET TO MOTOR WITH ARROWS ORIENTS AS SHOWN. SECURE COLLET TO MOTOR WITH 75 ± 5 IN-LBS OF TORQUE


 Medtronic		
SEE WORK INSTRUCTIONS FOR SPECIFICATIONS		
SCALE: 1:1	REV: 3 OF 4	REV: 01/11/11
TITLE		
ASST. FINAL, LEGEND GOLD TOUCH		
DESIGN NUMBER 2175-0003		REV P

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
DETAIL B
SCALE 3:2
ONCE ASSEMBLED:
 - ITEMS 19 AND 12 SHALL SWIVEL SMOOTHLY
 - ITEMS 11 AND 12 SHALL SWIVEL SMOOTHLY WITH DEFINITE STOPS AT

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INDEX	QTY.	PART NUMBER	DESCRIPTION
1	1	100-075-20	MOTOR SEAL 2
2	1	100-510-05	ROTOR SHAFT-3 VANE OFFSET, LEGEND GOLD, RULON J SEAL
3	1	100-510-08	ROTOR HOUSING, LEGEND GOLD MOTOR
4	1	100-510-17	FERRULE, EXHAUST, MOTOR A
5	1	100-510-18	DISTRIBUTOR, MOTOR A
6	3	100-510-31	VANE, LEGEND GOLD MOTOR
7	1	100-510-38	REAR BEARING HOUSING, LEGEND GOLD MOTOR
8	1	100-510-50	SWIVEL CONVERTIBLE, HOSE END, MOTOR A
9	1	100-510-54	FRONT SWIVEL UNION, LEGEND GOLD MOTOR
10	1	100-510-55	COMPRESSION SHIM, MOTOR A
11	1	100-510-56	NUT, COMPRESSION, LEGEND GOLD MOTOR
12	1	100-510-57	NOZZLE, FRONT, MOTOR A
13	1	100-510-63	SOUND TUBE SPACER, GOLD MOTOR
14	1	100-510-64	HOUSING, FRONT BEARING, LEGEND GOLD MOTOR
15	1	100-510-66	BARB, HP HOSE, GOLD, CRIMPED
16	1	100-550-1	CASE, MOTOR, LEGEND GOLD
17	1	MT10423B001	ASSEMBLY, DOUBLE D COLLET, PNEUMATIC
18	1	120-030-1	OIL, MPSS REX LUBRICANT
19	1	120-051-3	EXHAUST TUBE, OILER
20	1	120-051-17	1.0 LONG EXHAUST HOSE, LEGEND PLATINUM
21	1	MT08764B001	ASSEMBLY, HOSE, HIGH PRESSURE
22	1	120-051-23	SPRING, COMPRESSION, MODIFIED, .430 OD X .041 WIRE X 10" LONG
23	1	120-051-24	SEAL, LIP, OILER HOUSING
24	1	120-051-25	STEM, LUBRICATOR/DIFFUSER, LEGEND
25	1	120-051-26	SCREEN, LUBRICATOR/DIFFUSER, LEGEND
26	1	120-051-28	CONNECTOR SPACER, OILER HOUSING, LEGEND
27	1	M405204B001	EXHAUST HOSE, LEGEND
28	1	120-051-30	PRESSURE TUBE, OILER, CRIMPED
29	1	120-051-33	HOUSING, OILER
30	1	120-057	ASSEMBLY, CAPILLARY TUBE, LEGEND GOLD

31	1	900-188	BEARING, .250 ID X .500 OD
32	1	900-2-5	BEARING, .125 ID X .3125 OD
33	1	901-001	BALL, .078 DIA, 440C, SS
34	56	901-009	BALL, 3/64" (.0469) DIAMETER
35	1	902-013	PIN, SPIRAL, .049 X .500
36	1	903-012	O-RING, .003 CS, 0.789 ID (FLUOROELASTOMER)
37	1	903-016	O-RING, .031 X-SEC X .312 ID, 70 DURO SILICONE
38	1	903-017	O-RING, .032 X-SEC X .062 ID, 70 DURO SILICONE
39	1	903-030	O-RING, .079 CS, X .945 I.D.
40	1	903-037	O-RING, .032 X-SEC, .220 I.D., 70 DURO SILICONE
41	1	903-038	O-RING, .040 X-SEC, .226 I.D., 70 DURO SILICONE
42	1	903-039	O-RING, .040 X-SEC, .226 I.D., 70 DURO SILICONE
43	1	903-040	O-RING, .051 X-SEC, .547 I.D., 70 DURO SILICONE
44	1	903-30-622	O-RING, .030 CS X .622 ID
45	1	905-086	SPRING, .062 O.D. X .015 WIRE X .160 LONG
46	1	906-050	PLUG, 1/4" MALE BSPP
47	2	906-063	FERRULE, BRASS, .312 ID X .50 L X .225 HOLE, HPH
48	1	910-001	LOCTITE 680
49	3	910-004	WIRE, STAINLESS STEEL, .020 DIAMETER
50	1	910-014	LOCTITE 609
51	1	910-015	PRIMER "M"
52	1	910-017	SUPER LUBE SYNTHETIC GREASE W/PTFE
53	1	910-096	DC 732 RTV SILICONE ADHESIVE SEALANT (NOT SHOWN)
54	1	910-098	LOCTITE 545 THREAD SEALANT
55	1	910-120	SLEEVE, NOISE REDUCTION, LEGEND MOTORS
56	1	911-003	SET SCREW, 2-56 X .093
57	2	911-006	SCREW, SET, SOCKET HEAD, 2-56 X 1/16"
58	2	911-041	SET SCREW, SOCKET HEAD, 4-40 X 1/16"
59	3	912-036	SCREW, #4-40 X 3/16 SS SOCKET HEAD CAP
60	1	913-004	SCREW, PAN HEAD, 2-56 X .125

THIS DRAWING IS A COPY OF THE ORIGINAL DRAWING AND IS NOT TO BE USED FOR REPRODUCTION WITHOUT THE WRITTEN CONSENT OF THE COMPANY	SPECIFICATIONS (UNLESS OTHERWISE NOTED) DIMENSIONS IN INCHES UNLESS NOTED DECIMALS ALL SHOWN DIMENSIONS ALL FINISHES UNLESS NOTED FINISHES: IF NOT SPECIFIED SURFACES: UNLESS OTHERWISE NOTED: P.A. DIMENSIONAL TOLERANCES: FRACTIONS: .005 DECIMALS: .005 ANGLES: .005 HOLE DIA: .005 TAP DIA: .005 TAP DIA: .005	THIRD ANGLE PROJECTION DIMENSIONS AND TOLERANCES FOR ASME Y14.5 - 2009 ENGINEER: CHECKED: DATE: 1-10-2010 DATE: 1-10-2010	DESCRIPTION OF CHANGE TITLE ASSEMBLY, FINAL LEGEND GOLD MOTOR DRAWING NUMBER 2175-0001 SHEET P
---	--	---	--

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
[EDO] = ESSENTIAL DESIGN OUTPUT

○ = INSPECTION REQUIRED


- 1) ENSURE PARTS ARE FREE OF BURRS AND FIT/FUNCTION WELL. CLEAN ALL METAL COMPONENTS TO REMOVE ANY OIL, DIRT, GREASE OR DEBRIS BEFORE ASSEMBLY AS NEEDED.
- 2) ASSEMBLE PER WORK INSTRUCTIONS SPECIFIED IN ACCOMPANYING ROUTER.
- 3) APPLY MIDAS REX LUBRICANT (ITEM 14) TO O-RINGS PRIOR TO INSTALLATION.
- 4) APPLY PRIMER (ITEM 45) TO THREAD AND ALLOW TO DRY. APPLY ADHESIVE (ITEM 44) AND INSTALL.
- 5) APPLY THREAD SEALANT (ITEM 48) TO MALE THREADS OF SPACER (ITEM 22) AND INSTALL TO OILER HOUSING (ITEM 25).
- 6) DRILL $\varnothing .070 \times .046$ INTO CASE (ITEM 5). TAP #2-56 HOLE.
- 7) BRAID SHALL BE WITHIN 0.15" FROM HOSE END BEFORE CRIMPING
- 8) FERRULE SHALL BE PRESSED FLUSH AGAINST FLAT ON BARB OR PRESSURE TUBE.
- 9) COAT WITH SILICONE ADHESIVE (ITEM 47) TO PROTECT FROM WIRE ENDS AND SECURE HOSE.
- 10) LUBRICATE SEAL (ITEM 2) WITH LUBE (ITEM 46) PRIOR TO INSTALLING ROTORSHAFT (ITEM 3).
- 11) REFERENCE 600-023 FOR DEFINITION AND DELINIATION OF MIDAS TRUE TYPE FONT.
- 12) APPLY PRIMER (ITEM 45) TO THREAD AND ALLOW TO DRY. APPLY ADHESIVE (ITEM 41) AND INSTALL.
- 13) QUANTITY BASED ON TYPICAL DIMENSIONS. ADJUST FOR BEST PERFORMANCE.
- 14) LOCATE APPROXIMATELY AS SHOWN.

DETAIL D AND SERIAL NUMBER MARKING

SCALE 2:1

PRIMER (ITEM 45) THREADS OF COMPRESSION NUT. APPLY ADHESIVE (ITEM 44). ASSEMBLE COLLET TO MOTOR WITH ARROWS ORIENTS AS SHOWN. SECURE COLLET TO MOTOR WITH 75 ± 5 IN-LBS OF TORQUE

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
TRADE SECRET - DO NOT COPY
THIS MATERIAL IS CONFIDENTIAL AND PROPRIETARY, BEING
UNLAWFULLY OBTAINED BY MICROAIR INSTRUMENTS
FORT WORTH, TEXAS 76137. THE DISSEMINATION OF THE
INFORMATION CONTAINED HEREIN IS PROHIBITED BY FEDERAL
COPYRIGHT LAWS AND ALL GOVERNMENT PROSECUTIONS
AND MAY BE SUBJECT TO CRIMINAL PENALTIES.
DOCUMENT # MPSE-A-002 REV. 1

CROSS SECTION VIEW

NOTES:

1. SHIP UNIT WITH CERTIFICATE OF CONFORMANCE.
2. ETCH **Medtronic**
AT 0.25 OVERALL HEIGHT
CENTERED AS SHOWN
3. ETCH * CE 0086 SYMBOL *
AT 0.40 OVERALL HEIGHT
AS SHOWN
4. ETCH "SN *****" IN
0.06 HIGH HELVETICA
FONT AS SHOWN
5. CENTER WITH BLUE ANODIZE
PER SPECIFICATION SC602030
IN 0.20 HIGH SCALASANS
LINING FONT AS SHOWN,
TWO PLACES:

TRITON



TM
6. CENTER WITH BLUE ANODIZE PER SPECIFICATION SC602030
IN 0.05 HIGH SCALASANS LINING FONT AS SHOWN TWO PLACES:

TM
7. CENTER WITH BLACK ANODIZE PER SPECIFICATION SC602024
IN 0.08 HIGH HELVETICA FONT AS SHOWN TWO PLACES:

HIGH TORQUE HANDPIECE
8. CENTER WITH BLACK ANODIZE PER SPECIFICATION SC602024
IN 0.06 HIGH HELVETICA BOLD FONT AS SHOWN:

DO NOT IMMERGE
OR OIL
PN 720001
9. CENTER WITH BLACK ANODIZE PER SPECIFICATION SC602024
IN 0.15 HIGH HELVETICA BOLD FONT AS SHOWN

A. DO NOT IMMERGE
OR OIL
PN 720001
10. CLEAR (NATURAL) ANODIZE ALUMINUM BODY, EXCEPT TEXT
AND LOGO, PER SPECIFICATION SC602029
11. VENDOR LABEL M720001 PRINTED ON LABEL 170012ML (NOT SHOWN)
12. PLACE (FC 912008ML) (NOT SHOWN) IN BOX

VIEW A-A
UNFOLDED

VENDOR

MICROAIR INSTRUMENTS
CHARLOTTEVILLE, VA 22911

SPECIFICATIONS

TRITON HT HANDPIECE
PN 720001

FINAL RELEASE
2-2-01

ORIGINAL

REV	1	DATE	02/01/01	BY	DAVE	CHKD	DAVE	APPD	DAVE
SPECIFICATIONS (UNLESS OTHERWISE NOTED)									
ALL DIMENSIONS ARE IN INCHES									
BREAK ALL SHARP EDGES: .005 MAX									
ALL INSIDE CORNERS: .010 MAX									
FRESH V. MARK, QUANTITY MARKING									
TOLERANCES									
FRACTIONAL									
DECIMAL									
ANGLES									
SURF. FINISH									
MATERIAL: SEE VENDOR SPEC									
FINISHES: 2/2									
CREATION DATE: 4/14/00									
DO NOT SCALE FROM PRINT									
SHEET 1 OF 1									
SCALE: FULL									
DRAWING NUMBER									
720001									
REV									
F									

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Motors	Required Operating (Dynamic) Pressure	Nominal Operating (Dynamic) Pressure	Approximate Flow Rate Required	Gas Type
Pneumatic Motor Midas Rex Legend Gold	80 – 120 psi (551.6-827.4 kPa)	100 psi (689.5 ± 20 kPa)	12 cubic feet/min (339.6 liters/min)	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar (550-830 kPa)	6.9 bar (690 ± 20 kPa)	340 liters/min	
Pneumatic Motor Midas Rex Legend Gold Touch	80 – 120 psi (551.6-827.4 kPa)	100 psi (689.5 ± 20 kPa)	12 cubic feet/min (339.6 liters/min)	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar (550-830 kPa)	6.9 bar (690 ± 20 kPa)	340 liters/min	
Pneumatic Motor Midas Rex Legend Platinum	80 – 120 psi (551.6-827.4 kPa)	100 psi (689.5 ± 20 kPa)	12 cubic feet/min (339.6 liters/min)	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar (550-830 kPa)	6.9 bar (690 ± 20 kPa)	340 liters/min	
Pneumatic Motor Midas Rex MR7 Motor	80 – 120 psi (551.6-827.4 kPa)	100 psi (689.5 ± 20 kPa)	12 cubic feet/min (339.6 liters/min)	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar (550-830 kPa)	6.9 bar (690 ± 20 kPa)	340 liters/min	
Pneumatic Motor Midas Rex MR7 Touch Motor	80 – 120 psi (551.6-827.4 kPa)	100 psi (689.5 ± 20 kPa)	12 cubic feet/min (339.6 liters/min)	Nitrogen or Dry-Filtered Compressed Air
	5.5 – 8.3 bar (550-830 kPa)	6.9 bar (690 ± 20 kPa)	340 liters/min	
TRITON Motor High Torque Handpiece	120 – 200 psi (827.4-1379 kPa)	100 psi (689.5 ± 20 kPa)	12 cubic feet/min (339.6 liters/min)	Nitrogen or Dry-Filtered Compressed Air
	8.3 – 13.8 bar (830-1380 kPa)	7 bar (700 ± 20 kPa)	340 liters/min	

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Performance Specification:

Motors	Dimensional:	Performance			Duty cycle	Functional
		Speed	Torque	The maximum noise level		
Pneumatic Motor Midas Rex Legend Gold	Weight: 138 ± 10 g Height: 4.0 inches (10.16 ± 1 cm) Width: 0.73 inches (1.85 ± 0.1 cm) Thickness: 0.73 inches (1.85 ± 0.1 cm)	76 000-82 000 RPM @120psi (827.4 kPa)	4.7oz-9.6oz (0.133 – 0.27 kg) @120 psi (827.4 kPa)	77 dBA measured 18" behind the motor (0.46 m)	Designed for intermittent rotation: approximately 1 to 3 minutes of continuous use followed by non-use periods of 1 to 20 minutes.	The maximum allowed nose flow is 1.5 liters/minute.
Pneumatic Motor Midas Rex Legend Gold Touch	Weight: 138 ± 10 g Height: 4.0 inches (10.16 ± 1 cm) Width: 0.73 inches (1.85 ± 0.1 cm) Thickness: 0.73	76 000-82 000 RPM @120psi (827.4 kPa)	4.7oz-9.6oz (0.133 – 0.27 kg) @120 psi (827.4 kPa)	79 dBA measured 18" behind the motor (0.46 m)	Designed for intermittent rotation: approximately 1 to 3 minutes of continuous use followed by non-use	The maximum allowed nose flow is 1.5 liters/minute.

		inches (1.85 ± 0.1 cm)				periods of 1 to 20 minutes.	
	Pneumatic Motor Midas Rex Legend Platinum	Weight: 138 ± 10 g Height: 4.0 inches (10.16 ± 1 cm) Width: 0.73 inches (1.85 ± 0.1 cm) Thickness: 0.73 inches (1.85 ± 0.1 cm)	76 000-82 000 RPM @120psi (827.4 kPa)	4.7oz-9.6oz (0.133 – 0.27 kg) @120 psi (827.4 kPa)	77 dBA measured 18" behind the motor (0.46 m)	Designed for intermitted rotation: approximately 1 to 3 minutes of continuous use followed by non-use periods of 1 to 20 minutes.	The maximum allowed nose flow is 1.5 liters/minute.
	Midas Rex MR7 Pneumatic Motor	Weight: 138 ± 10 g Height: 4.0 inches (10.16 ± 1 cm) Width: 0.73 inches (1.85 ± 0.1 cm) Thickness: 0.73 inches (1.85 ± 0.1 cm)	66 000 – 82 000 RPM @120psi (827.4 kPa)	4.7oz-9.6oz (0.133 – 0.27 kg) @120 psi (827.4 kPa)	82.4 dBA @ 80-120 psi (827.4 kPa)	Designed for intermitted rotation: approximately 1 to 3 minutes of continuous use followed by non-use periods of 1 to 20 minutes.	The maximum allowed nose flow is 1.5 liters/minute.
	Midas Rex MR7 Touch Pneumatic Motor	Weight: 138 ± 10 g Height: 4.0 inches (10.16 ± 1 cm) Width: 0.73 inches (1.85 ± 0.1 cm) Thickness: 0.73 inches (1.85 ± 0.1 cm)	66 000 – 82 000 RPM @120psi (827.4 kPa)	4.7oz-9.6oz (0.133 – 0.27 kg) @120 psi (827.4 kPa)	82.4 dBA @ 80-120 psi (827.4 kPa)	Designed for intermitted rotation: approximately 1 to 3 minutes of continuous use followed by non-use periods of 1 to 20 minutes.	The maximum allowed nose flow is 1.5 liters/minute.
	TRITON Motor High Torque Handpiece	Weight: 11b 5.6 oz. (606.6 gm) Height: 5.72 in (14.53 cm) Width: 3.72 in (9.45 cm) Thickness: 1.08 in (2.74 cm)	1 600 RPM nominal output speed	125 in/oz (8 in/lb/ 0.88 H-m) at 300 RPM minimum @ 100 PSI (689.5 kPa)	82.4 dBA @ 80-120 psi (827.4 kPa)	Designed for intermitted rotation: approximately 1 to 3 minutes of continuous use followed by non-use periods of 1 to 20 minutes.	Trigger: Rotate for safe position Instrument to be cannulated to accept passage of a Dia 3.2 mm (0.125")
The maximum vibration allowed is 2.0 G The maximum temperature measured after five minutes of motor free running should be 105°F (40.55°C)							
	Принадлежности:		Размеры	Вес	Accessories:	Dimensions	Weight
	Картриджи смазочные Midas Rex Lubricant\Defuser, модель PA100-A		4,5" x 4,5" x 2,5" (11,4±0,5 x 11,4±0,5 x 6,4±0,2 cm)	52,5 ±5 г	Midas Rex Legend Lubricant/Diffuser Pack, PA100-A	4.5" x 4.5" x 2.5" (11.4±0.5 x 11.4±0.5 x 6.4±0.2 cm)	52.5 ±5 g
	Картриджи смазочные Midas Rex Lubricant\Defuser, модель PA700		8,25" x 2,68" x 3,88" (21±0,8 x 6,8±0,2 x 9,9±0,5 cm)	10 ±2 г	Midas Rex MR7 Lubricant/ Diffuser Pack, PA700	8.25" x 2.68" x 3.88" (21±0.8 x 6.8±0.2 x 9.9±0.5 cm)	10 ±2 g
	Ключ для мотора Midas Rex, модель PA130		5,0" x 0,75" x 0,120" (12,7±0,5 x 1,9±0,1 x 0,3±0,02 cm)	54 ±5 г			

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



TRITON Handpiece

The Handpiece it is placed in a folded carton tray which is then inserted into a poly bag. The bag is then placed inside a labeled box. The Important Information Brochure is placed in box.

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
Pneumatic Control Units

The control unit is placed in a folded carton tray which is then inserted into a poly bag. A label is attached to the outside of the bag. The bag is then placed inside a labeled box.


Regulator

The regulatory is placed in a folded tray and then inserted into a labeled box.

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Lubricants/Diffusers

Lubricants/Diffusers are placed in the box with the caps at the top of the box – 4 per box. Each PA700 individually placed in poly bag. The box is then closed and labeled appropriately.


Cleaning Brushes

Cleaning Brushes are placed in poly tube. The tube is then closed and labeled appropriately.

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


Information in Russian sticker

The final packaging for the Russian market is accompanied with additional sticker. Devices supplied in packaging marked with the following possible information:

Other accessories

Other accessories are inserted into a poly bag which is then heat sealed then labeled.

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1. Name of medical device in Russian language;
2. Name of the manufacturer;
3. Country of origin;
4. Code UPN (unique product number);
5. Code CFN (customer facing number);
6. Date and number of Registration Certificate;
7. Sign of conformity in Russia ;
8. Warning and informing symbols and inscriptions (if applicable).

Sample of the sticker

		
	Manufacturer:	Code UPN
	Country of origin:	Code CFN
	Manufacturer's authorized representative::	
	Registration Certificate № _____ от _____	
	Date of manuf./Lot/U D:	See IFU
	See on pack	See packaging symbols

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	Catalogue number
	Lot number
	Manufacturer
	Non sterile
	For single use only!
	Authorized representative in the European Community
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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running, or when the motor or attachment is in an overheated state.

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DO NOT use Legend Pneumatic Control Unit finger/foot control in the “ ” position to operate a Legend Gold or Legend Platinum

Motor. The control is used in the “ ” position only for the Legend Gold Touch Motor.

DO NOT install or remove an attachment and dissecting tool from the Legend Gold Touch unless the safety slide is in the “safe” position.

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- Do not use the Midas Rex MR7 System before proper cleaning and sterilization.

- Do not operate the Midas Rex MR7 System in the presence of Magnetic Resonance Imaging devices.

- Do not use damaged, faulty, or modified Midas Rex MR7 System components. Inspect the Midas Rex MR7 System for damage prior to each use:

- Check the motor's exhaust hose for cracks or tears.

- Visually inspect attachments and tools. Do not use bent or damaged tools.

- Install attachment and dissecting tool, then briefly run motor.

- Check motor for overheating and leaking lubricant.

- * Check attachment for overheating.

- * Check dissecting tool for I ail.

- Do not operate the Midas Rex MR7 System without eye protection.

- Motors and attachments which fail due to extended use may allow a component to detach and fall from the motor or attachment, and may cause patient injury.

- Heavy side loads and/or long operating periods may cause the device to overheat. If overheating occurs:

- Never place an overheated motor on the patient or draping during surgery.

- Discontinue use and rest the motor by using intermittently, or wrap the motor/attachment interface with a moist sterile towel.

- If the motor is passed of, the receiver should grasp the motor by the proximal end close to the motor hose.

- To avoid injury to the patient or user, do not place the handpiece on the patient or in an unsecured location, when not in use.

- Midas Rex MR7 motors should only be operated when the attachment is in the position.

If a dissecting tool package is opened, but the tool is not used or contaminated, the

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		1B. Flush the nose of the TRITON HT Handpiece, Reciprocating Saw Attachment and Compact Micro-Reciprocating Saw Handpiece with a water-pik.	
		2. Rinse instrumentation under running water to remove all traces of de ter gent solution. If possible, use dis tilled water for the final rinse.	
		3. Disconnect the TRITON Hose from the TRI TON HT Handpiece or Compact Micro-Saw Handpiece. Dry all components with a lint-free towel.	
		Remove TRITON™ Hose from TRITON™ Handpiece. DO NOT attempt to remove hose from the MIDAS REX® Motor or MEDNEXT® Drill.	
		Disinfection	No particular requirements
		Sterilization (Temperatures are minimum required, times are minimum required)	Remove TRITON™ Hose from TRITON™ Handpiece. DO NOT attempt to remove hose from the MIDAS REX® Motor or MEDNEXT® Drill.
		2. Place cleaned MIDAS REX® Motor or MEDNEXT® Drill and TRI TON™ Hose in bottom of instrument case.	
		3. Place TRITON™ Handpieces, attachments, and accessories in appropriate depressions of insert trays. Set tray in case.	
		NOTE: Midas Rex or Mednext attachments can be placed in their respective insert trays and positioned in the TRITON instrument case.	
		4. Place lid on instrument case and sterilize.	
Steam Sterilization:			
Cycle:	Gravity	Pre-vac	
Temperature:	132°C/270°F	132°C/270°F	
Time:	for 30 minutes.	for 4 minutes.	

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Please refer to the preventive maintenance manual for the specific steps necessary to maintain Midas Rex, Legend or TRITON System.

MANUFACTURER'S RECOMMENDATION

It is strongly recommended that surgeons, operating room personnel, and service personnel become thoroughly familiar with all Midas Rex, Legend or TRITON Systems prior to use in surgery. To assist in this effort, Medtronic Powered Surgical Solutions offers several Hands-On Workshops. Contact your Medtronic Neurologic Technologies Sales Representative or your Medtronic Regional Distributor for schedules.

TRITON SYSTEM TROUBLESHOOTING

When TRITON equipment requires servicing or repair, it is strongly urged that the equipment be returned to MEDTRONIC MIDAS REX, Fort Worth, TX, for quality assured service by factory-trained personnel who will utilize genuine parts as required.

LEGEND and MR7 SYSTEM TROUBLESHOOTING

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Motor not Running or Low on Power:

Possible Cause	Solution
Hoses not properly connected. Operating pressure inadequate.	Make sure all connections are secure. Check gas supply pressure gauge. Increase pressure according to compressed gas requirements, if necessary.
Attachment not properly installed and locked onto the motor.	Remove and re-install attachment and tool to ensure proper installation and locking of attachment onto motor.
Foot pedal on pneumatic control unit not functioning properly.	Check for obstructions under the foot pedal. If foot pedal continues to fail, return the pneumatic control unit to MPSS to be refurbished.
Motor stalls.	Manually spin the dissecting tool then activate the motor. If the motor continues to stall, return it to MPSS to be refurbished.

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	site during bone dissection.
Attachment Will not Properly Seat on the Motor:	
Possible Cause	Solution
Motor collet flats are not aligned.	Use the Legend motor wrench to rotate the flat closest to the motor case until its marker is aligned with the marker on the flat farthest away from the motor case.
Tool is Difficult to Remove from Attachment:	
Possible Cause	Solution
Aging of attachment. Improper cleaning. Use of reprocessed tools. Use of an unauthorized refurbisher.	Return to MPSS to be refurbished, or purchase new equipment.
16-MF Contra-Angle Attachment is Overheating	
Possible Cause	Solution
The Contra-Angle attachment operates by a set of internal gears to engage the drive shaft. It is normal for some heat to be generated approximately 2cm from the distal end of the attachment and at the right of the angle head.	Verify pressure setting of 80 psi (5.5 bar). If heat continues or is excessive, return the attachment to MPSS to be refurbished or purchase new equipment.
Perforator is Running too Slow	
Possible Cause	Solution
Pressure set incorrectly.	Check the pressure setting at the foot control.
Dissecting Tool Flails	
Possible Cause	Solution
A non-Legend dissecting tool is being used. Worn attachment or tube bearings.	Replace with a Legend dissecting tool. Try another attachment or tube to isolate the location of the problem. If the attachment is failing, return it to MPSS. If the tube is failing, dispose of it and use a new tube.
Attachment/tube and tool are not compatible.	Match color code on the dissecting tool packaging to the color code on the attachment/tube.
Motor is damaged. Dissecting tool's size and geometry may contribute to flailing at certain speeds.	Return motor to MPSS to be refurbished. Adjust the speed by changing the pressure setting or foot/finger control. Do not use if flailing persists. Change dissecting tools.
Dissecting Tool Vibrates Excessively	

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PROBLEM	POSSIBLE CAUSE	RECOMMENDATIONS
Midas Rex LEGEND MOTORS		
Motor does not run or is low on power.	Hoses not properly connected.	Check all connections from compressed gas source to motor.
	Operating pressure inadequate.	Check gas supply to Pneumatic Control Unit at Supply Pressure gauge. Pressure may need to be increased according to compressed gas requirements.
	Attachment not properly installed and locked onto the motor (possible smoking condition).	Remove and re-install attachment and dissecting tool to ensure proper installation and locking of attachment onto motor.
	Motor sterilized with single use, disposable Lubricant/Diffuser cartridge attached to the motors lubricant/diffuser housing.	DO NOT use. Return motor to MPSS to be refurbished.
	Foot Pedal on Pneumatic Control Unit not properly functioning.	Check for obstruction under foot pedal. If foot pedal continues not to

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		function properly, return the Pneumatic Control Unit to MPSS to be refurbished.
	Motor stall or will not spin.	Manually spin the dissecting tool then activate the motor. If motor continues to stall, then return motor to MPSS to be refurbished.
Legend Gold Touch will not run after checking possible causes listed previously for all Legend motors.	Finger/Foot control on the Pneumatic Control Unit not activated.	Push finger/foot control on Pneumatic Control Unit to the  position.
	Safety slide in  position.	Move safety slide to  position.
Motor continues to run.	Pneumatic Control Unit is not properly functioning.	Depressurize system. Return Pneumatic Control Unit to MPSS to be refurbished.
	Finger control is not properly functioning	Return motor to MPSS to be refurbished.
Motor is missing a component or a component has detached from the motor.	Motor failed.	DO NOT use. Return motor to MPSS to be refurbished.
Motor rotates, but an abnormal noise is heard.	Lubrication inadequate.	Check for proper installation of the Legend Lubricant/Diff user. If problem persists, replace Legend Lubricant/Diffuser. If problem continues, return motor to MPSS to be refurbished.
	Motor's exhaust hose damaged or cut, or internal pressure hose detached.	DO NOT use. Depressurize system. Return motor to MPSS to be refurbished.
	Bearings of motor or attachment worn.	Change attachment to isolate location of problem to the attachment or motor. Return motor and/or attachment to MPSS to be refurbished.

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defect in materials or workmanship within a period of 90 days from the date of shipment of a new Accessory Component or serviced or refurbished Accessory Component to the purchaser, Medtronic will at its option: (a) issue a credit to the purchaser of the replacement Motor Component or Accessory Component equal to the Purchase Price, as defined in Subsection A(3), against the purchase of any same Motor Component or Accessory Component requested as its replacement, or, (b) provide a functionally comparable replacement Motor Component or Accessory Component at no charge. (3) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original or current functionally comparable, or replacement Motor Component or Accessory Component.

B. To qualify for this Limited Warranty, these conditions must be met:

- (1) The Legend Pneumatic High Speed System must be used in accordance with the labeling and Instruction Manual provided with the products.
- (2) Only Legend Pneumatic High Speed System dissecting tools that are manufactured, authorized, and quality controlled by Medtronic may be used.
- (3) No reprocessed dissecting tools may be used.
- (4) The Legend Pneumatic High Speed System may not be damaged by accident, abuse, misuse or tampering.
- (5) The Legend Pneumatic High Speed System may not be modified without the written permission of Medtronic.

Triton

C. This Limited Warranty is limited to its express terms. In particular:

- (1) Except as expressly provided by this Limited Warranty, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE LEGEND PNEUMATIC HIGH SPEED SYSTEM TO FUNCTION WITHIN NORMAL TOLERANCES WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, NEGLIGENCE, STRICT LIABILITY, OR OTHER TORT OR OTHERWISE.
- (2) This Limited Warranty is made only to the purchaser. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM OR OTHERWISE.
- (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Limited Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid. This Limited Warranty gives the purchaser specific legal rights. The purchaser may also have other rights that vary from state to state.

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damage caused by the use of non-genuine MIDAS REX, MEDNEXT or TRITON™ manufactured cutting tools. Peak performance, reliability and maximum service life from your MIDAS REX, MEDNEXT or TRITON™ instrumentation may be assured by using only those products that are manufactured by and sold through Medtronic Midas Rex, Fort Worth, Texas.

Midas Rex MR7

A. This Limited Warranty provides the following assurance to the purchaser of a Medtronic Midas Rex MR7 Pneumatic High Speed System. This Limited Warranty is extended only to the buyer purchasing the MR7 System directly from Medtronic or from its affiliate or its authorized distributor or representative. The Midas Rex MR7 Pneumatic High Speed System includes the motor, foot control, instrumentation cases and trays (hereafter referred to as System Components), straight and angled motor attachments (hereinafter referred to as "Attachments"), telescoping tubes (hereinafter referred to as Semi-reusable Components) and dissecting tools and other accessories not listed above and jointly referred to as MR7 Pneumatic High Speed System, unless specifically noted.

(1) Should a System Component fail to function to Medtronic's published specifications during the term of this Limited Warranty (one year from the date of sale of a new System Component or 90 days from the date of sale of a refurbished or used System Component), Medtronic will either repair or replace the Motor Component or any portion thereof.

(2) Should an Attachment fail to function to Medtronic's published specifications during the term of this Limited Warranty (90 days from the date of sale of a new Attachment), Medtronic will either repair or replace the Attachment or any portion thereof.

(3) Should a Semi-reusable Component fail to function to Medtronic's published specifications during the term of this Limited Warranty (30 days from the date of sale of a new Semi-reusable Component), Medtronic will replace the Semi-reusable Component or any portion thereof.

(4) Should a Single Use Component fail to function to Medtronic's published specifications prior to its "use by" date Medtronic will replace the Single Use Component. B. To qualify for this Limited Warranty, these conditions must be met:

(1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.

(2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

(3) Medtronic must be notified in writing within thirty (30) days following discovery of a defect.

(4) The Product must be returned to Medtronic within thirty (30) days of Medtronic receiving notice as provided for in (3) above.

(5) Upon examination of the Product by Medtronic, Medtronic shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and

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services, if applicable, have been performed on the Product

C. This Limited Warranty is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the MR7 System, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this Limited Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

Transmissible Spongiform Encephalopathies (TSE) Return Policy

Medtronic Powered Surgical Solutions (MPSS) will not authorize or accept the return of MPSS products that directly contact patients or is contaminated with a patient's body fluids who is suspected or confirmed with a Transmissible Spongiform Encephalopathies/Creutzfeldt-Jakob Disease (TSE/CJD) diagnosis. Furthermore, MPSS recommends that all MPSS products used on a patient confirmed with a TSE diagnosis be incinerated. Contact your Sales Representative for replacement of product incinerated under this policy or for temporary equipment while original equipment is quarantined.

Contact MPSS Regulatory Affairs Department for additional information regarding TSE contamination. MPSS dissecting tools used on a patient suspected of a TSE/CJD diagnosis must be incinerated.

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