Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuangxin Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.,
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 (0)24 6258 9911
Fax: +86 (0)24 6258 9922
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseaa.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseaa.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 5250-5607, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquin Floriano, nº. 413, 20º andar/floor. Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskopi Sverige AB
Storsatragrand 14,
12724 Skärholmens, Sweden
Postal address: P.O. Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Mälahdenkatu 23a, 00210 Helsinki, Finland
Phone: +35-896824774, Fax: +35-8968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopleer.be

KARL STORZ Endoscopia France S. A.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colanin, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Gounari 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantiou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maitepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

OOO KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: info-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
Solidiere - Beirut Souks
Block M, 3rd Floor
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061
Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Bokeykhan Ulica 8/2, VP - 1
010000 Astana, Respublika Kasachstan
Phone: +7 7172 57-52 16 / 57 28 49 / 57 09 34
Fax: +7 7172 43 96 96
E-Mail: info@karlstorz.kz

KARL STORZ Endoscopy - Gulf & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O.Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.ksegn@karlstorz-roe.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 26812986
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F 10-0, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuangxing Mansion,

No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911
Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company,
Room 1119-20, Dongshan Plaza
69 Xianle Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
06-07 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 65325548, Fax: +65 65323832
E-Mail: infoSGP@karlstorz.com.sg
E-Mail: servicesSGP@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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**Important information for users
of KARL STORZ units and
instruments**

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This manual contains information for personnel authorized to prepare instruments for proctology. The instruments and accessories manufactured by KARL STORZ are, as a rule, sensitive precision instruments which require particular care during preparation. To prevent damage during preparation, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ endoscopic instruments.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply, design, features and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations and descriptions in this manual.

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**Important information for users
of KARL STORZ units and**

instruments U2

About these instructions	2
Using this manual	2
Symbols and hazard information	3
Safety instructions	4
General instructions for cleaning and sterilization	5

Rectoscopy 9

Rectoscope tube	9
Magnifier	11
Illumination rod	12
General instructions for forceps and scissors	13
CLICKline instruments with and without connector pin for unipolar coagulation	15
Take-apart instruments with U-spring handle	18
General instructions for coagulation instruments	20
Suction-coagulation cannulas	22
Adaptor for coagulation and dissecting electrodes	22
Suction and irrigation tubes with two-way stopcock	23
Instruments for polypectomy	22

Proctoscopes 25

Pediatric and adult proctoscopes	25
Anal retractor and speculum	26
Suction ligature instrument	27
Video operating rectoscope	28

Maintenance 29

Operating/storage conditions	29
Disposal	29
Directive compliance	29
Limitation of liability	29
Manufacturer's warranty	30
Important information	30
Servicing and repair	31
Repair and replacement program	31

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Using this manual

This manual is not an instruction manual, but instead provides general information regarding cleaning, disinfection, care and sterilization of proctology instrument sets from KARL STORZ.

This manual presumes knowledge of the "Cleaning, sterilization and care of KARL STORZ instruments" manual, in which the various cleaning and sterilization methods are described in detail. If you do not yet have this manual, it can be requested without charge from KARL STORZ. Instruments whose cleaning, dismantling and reassembly do not require any special knowledge are not listed separately in these instructions.

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Safety information

WARNING: The listed cleaning, disinfection and sterilization methods alone do not guarantee disinfection or sterility. This can only be achieved if staff follows recognized and validated preparation methods.

WARNING: New instruments are delivered unsterilized and must therefore be thoroughly cleaned, disinfected and/ or sterilized prior to initial use and each subsequent reuse.

WARNING: For all cleaning of contaminated endoscopes and instruments, the Accident Prevention Regulations (UVV) must be observed (e.g. wearing of protective gloves and safety goggles).

WARNING: Remember that sterilization is only then successful if

- the instruments have been freed of all organic materials and cleaning residues, as only clean surfaces can be reached by the sterilization.
- work proceeds according to recognized and validated methods.
- the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.

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CAUTION: Use of preparation methods other than those listed here may have a negative effect on the function and service life. Any damage resulting from this is not covered by the guarantee.

CAUTION: The following information on cleaning and sterilization only applies to the thermostable instrument set described in these instructions. For telescopes and fiberscopes, special rules apply (see instructions for telescopes and fiberscopes).

CAUTION: For the production and application of solutions, precisely observe the information from the manufacturer regarding mixing ratios and immersion times. A maximum immersion time of 60 minutes (including in water) must not be exceeded. Prolonged immersion can lead to damage. After cleaning it is recommended to rinse the instruments with completely demineralized water.

CAUTION: Constant changing between various preparation methods is especially hard on the materials and must therefore be avoided. A one-off change of method, due to procurement of new cleaning and sterilization equipment for example, is however possible.

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
General instructions on cleaning and sterilization

The medical device cycle shown below illustrates the sequence of steps for advisable instrument preparation:

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Used and contaminated instruments must be prepared immediately. Jointed instruments and scissors must be treated in the open state. For manual cleaning and disinfection, it is essential to follow the manufacturer's instructions for concentration and exposure time. Never use metal brushes or metal sponges for manual cleaning.

All scissors, coagulation instruments and micro-instruments which cannot be dismantled, should undergo preliminary ultrasound cleaning (exposure time of 3 - 5 minutes). Hollow instruments (suction devices etc.) should be held at an angle in the ultrasound bath so that the hollow spaces can fill. Position large instruments in the bath so that no acoustic shadow is formed. The following items, amongst others, are not suitable for cleaning in an ultrasound bath: all rigid and flexible telescopes, electrical devices, cables, prisms, rubber articles, magnetic discs and optical waveguides.

After cleaning, rinse thoroughly with completely demineralized water and then dry immediately. Rinse out hollow instruments with a cleaning pistol and then dry with compressed air. If necessary, give instruments such as drills and burs subsequent treatment with a cleaning brush.

If possible, cleaning disinfection should be carried out by machine, followed by steam sterilization at 134°C. For cleaning and disinfection, please observe the general information contained in the appendix and also the list of products approved by KARL STORZ for the particular purpose.

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All instruments should be checked for correct functioning before sterilization at the latest. This functional check should include the following aspects in particular:

- Smooth operation of the joints
- General signs of wear
- Exact closing of forceps, clamps and punches etc.
- Secure closing of valves and seals
- Sharpness of cutting instruments
- Condition of telescopes
- Condition of the insulating sheaths of cables, including cold light cables
- Broken fibres in cables

Damaged instruments must be sorted out and sent for expert, qualified repair. For this purpose we recommend the KARL STORZ Repair Service, which ensures consistent quality and guaranteed safety of your instruments and units.

Instruments with closure or joint must be treated with steam-sterilizable protective agents (e.g. instrument oil Art. no. 27656 B) before sterilization. As there are various types of cleaning, disinfection and sterilization, it is essential to ensure that recognized and validated (tested) preparation methods are used. Refer to the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ instruments."

CAUTION: The instructions for cleaning and sterilization in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ instruments" must be observed. There the methods for cleaning, disinfection and sterilization are explained in detail.

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
Rectoscope tube


Disassembly:

Take off the Illuminator head unit (turn the locking lever counterclockwise).

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
Take both adapters off the light port.

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Cleaning

NOTE: For heavy soiling, we recommend using an ultrasound bath for up to five minutes.

CAUTION: Do not clean the cap with glass window or the light carrier in an ultrasound bath.

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Manual cleaning

Remove tissue remnants with a cleaning brush. Brush and rinse out sheath tubes on the inside. Blow dry with sterile compressed air.

Machine cleaning


Make certain all parts are positioned securely to prevent damage during cleaning. The inside of the lumen must be rinsed out.

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
The glass surface of the light carrier must be carefully dried with alcohol after cleaning. Residues of disinfectants and cleaning agents on the glass surfaces may burn on when the optical waveguide is connected, and this inhibits light transmission considerably.

Firmly adhering coatings on the glass surfaces can be removed with special cleaning paste (Art. no. 27661).

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Assembly

Push the illuminator head unit onto the tube so that the red markings are opposite each other. Close the locking lever by turning clockwise.

Push the cap onto the illuminator ring.

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Changing the window of the cap

Take off the screw ring (turn counterclockwise).

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Take the window and seal ring off the cap.

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
Insert the new window and, if necessary, a new sealing ring.

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Magnifier

The magnifier is mounted on the illuminator head unit instead of the cap.

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
An instrument can be introduced into the working channel through the sealing cap ①. To introduce instruments with a larger diameter, the insert is removed.

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Cleaning/disinfection:


Only wipe-down disinfection is permitted for the magnifier.

Pull the sealing cap ① off the insert and, if necessary, replace it with a new one.

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
To replace the O-ring, undo the thumb screw, take out the insert and remove the O-ring.

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
Illumination rod

Before cleaning, unscrew and remove both the adaptors. The illumination rod can be autoclaved at up to 134 °C, gas or plasma sterilized or disinfected in a solution.

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The glass of surface the light carrier must be carefully dried with alcohol after cleaning. Residues of disinfectants and cleaning agents on the glass surfaces may burn on when the optical waveguide is connected, and this inhibits light transmission considerably. Firmly adhering coatings on the glass surfaces can be removed with special cleaning paste (Art. no. 27661).

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CAUTION: Do not clean the illumination rod in an ultrasound bath.

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General information for forceps and scissors

Disassembly and preparation for cleaning

Preparation should be started as soon as possible after using the instruments, since incrustations may cause irreparable damage. Preparation must begin no later than 8 hours after use.

To prevent damage during transportation and storage, the instruments must be carefully placed in suitable holders.

Cleaning

Clean surfaces, inside cavities and jaws thoroughly with a suitable brush, sponge or cleaning gun. Pay particular attention to the jaws when cleaning, since these are precision mechanical parts.

All scissors and forceps which cannot be completely dismantled must be cleaned ultrasonically.

If an irrigation port is available then the instrument must be rinsed out using this. For this purpose either connect the instrument to the water supply via the LUER lock connector or rinse out the irrigation channel using the cleaning gun (for at least 1 minute).

After cleaning, rinse thoroughly and dry with sterile compressed air. Instrument channels are especially easy to rinse with water using a cleaning gun and accessories (Art. no. 27660) and dried with sterile compressed air.

Manual cleaning

① NOTE: For heavy soiling, we recommend using an ultrasound bath for up to five minutes.

Remove sticky or caked tissue remnants on the jaws with a cleaning brush.

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CLICKline Instruments

CLICKline instruments can be dismantled very easily and may therefore be cleaned especially simply and thoroughly. Due to its modular design, the user may choose at any time among various handles. Furthermore, in case of damage to heavily used parts such as the jaws, only the defective element must be replaced instead of the complete instrument.

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CLICK'line instruments with and without connector pin for unipolar coagulation

Disassembly

Open handle completely. If necessary, unlock the locking mechanism on the handle.

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Press and hold the locking button ⑤ (do not turn).
The handle is now opened completely in order to pull out the sheath tube ②.

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
To remove the internals, hold the jaws ① tightly and open the bayonet mount ⑦ by rotating the sheath tube ② by a ¼ turn. Pull out the internals.

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
CAUTION: Please note that in instruments in the 2 and 3 mm sizes and in curved instruments the internals cannot be removed from the sheath tube.

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Manual cleaning

Remove tissue remnants with a cleaning brush (Cat. no. 27650 G). Brush and rinse out sheath tubes on the inside. Blow dry with sterile compressed air.

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

Machine cleaning/Disinfection

Make certain all parts are positioned securely to prevent damage during cleaning. The sheath tube must be rinsed out on the inside. The LUER-lock irrigation connection ④ permits easy attachment to a cleaning machine.

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
Assembly

Check all parts for damage. Replace if necessary. Check jaws and joints especially carefully, and additionally the insulation on instruments for HF use. Insert the internals into the sheath tube ② from the front. Lock the internals by rotating a ¼ turn in the bayonet mount on the sheath tube ②.

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Insert the sheath with the jaws closed as far as it will go into the completely opened handle.

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Then lock the sheath tube by closing the handle completely.

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
Place a drop of instrument oil (Art. no. 27656 B) in the joints and open and close the jaws several times.

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
In the case of handle 33140, the thumb grip can be removed. To do this, pull the thumb grip toward distal and remove by screwing in a clockwise direction.

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
When using handles 33134/33135, lubricate the joints and moving parts with a drop of instrument oil (Art. no. 27656 B). Test for proper operation.

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
**Take-apart instruments with
U-spring handle
Disassembly**

Press the guide wire ① downward and then detach from the grip.

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Cleaning

Clean the surfaces and jaws with a sponge or brush.

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

Assembly

Reassemble the instrument in reverse order.

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The coagulation surfaces of coagulation instruments should be wiped off with a moist surgical pad between uses during surgery. Incrustations may be loosened by briefly immersing in a 3% H₂O₂ solution. Before reusing, rinse with sterile water and wipe off again.

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NOTE: The HF cable and electrode connectors must be clean and dry. Moisture and dirt or residue from disinfectants may lead to malfunctions during operation.

Assembly and care

Incidents during high-frequency surgery which have resulted in damage to the instrument set or even injuries to the user or patient, may be traced back, under closer examination, to incorrect use or defective equipment. A frequent source of errors here proved to be damaged cables and/or electrodes. A further cause may lie in the use of unsuitable accessories from other manufacturers which are not compatible with KARL STORZ units. Safe and proper use of KARL STORZ instruments is only possible using KARL STORZ high-frequency cables and electrodes.

When testing for proper operation, check the insulation material especially for damage such as cracks and chipped material.

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Cable connections

WARNING: The use of damaged HF cables may lead to very serious dangers. Check the cable before every use for visible damage. Damaged cables must never be reused.

HF cables are subject to a certain degree of wear, particularly with improper handling and maintenance. For example, damage to the cable terminal or cable insulation may lead to HF voltage flashovers on the working element and therefore to the metallically conductive parts of the instrument. The parts of the instrument at HF potential then become dangerous to the patient and surgeon during patient application.

Careful inspection of the HF cable used is therefore absolutely necessary before every application.

Particular attention must be paid to

- burn marks on the cable terminal
- defects at the connections
- breaks in the cable or cable insulation
- missing pin for closing the second hole on the cable terminal.

If damage is found or suspected on HF cables, it is essential to discard them.

Cleaning high frequency cables

Clean the high frequency cable with a brush or sponge, or thermally disinfect in a cleaning and disinfection device. All KARL STORZ high frequency cables may be autoclaved up to 134°C, gas sterilized or immersed in disinfection solution. After immersing the cable in solution carefully dry it out in order to avoid any malfunctions due to residual moisture at the mating contacts. After cleaning, steam sterilization should preferably be performed at 134°C by the fractionated pre-vacuum method.

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Coagulation suction tubes Assembly

Check the seal ring ① for damage and, if necessary, replace with a new one.

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
Adaptor for coagulation and dissecting electrodes

To use 5 mm-sized coagulation or dissecting electrodes on handles intended for 10 mm suction/irrigation tubes, an adaptor is required.

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
Suction and irrigation tubes with two-way stopcock Disassembly

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Assembly

Lightly grease the contact surfaces on the cock plug ②. Insert the cock plug ② and fix with the spring cap ①. Tighten slightly. Test for proper operation.

NOTE: Due to hygiene-related reasons, the greased stopcock must be dismantled again for steam sterilization.

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


Unscrew and remove the HF adaptor from the socket ② and pull back the finger grip. Pull out the pull rod ① toward proximal.

Remove the remaining pieces of the loop from the pull rod.

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
Inserting the snare wire

Insert the pull rod ① into the guide tube from distal. Push the finger grip forward until the pull rod is located in its seat. Screw the HF adaptor into the socket ②.

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
Push the handle section forward as far as it will go. Insert the ends of the loop through the holes on the loop head from both sides so that the loop lies in the groove, and pull back the handle section 1. The loop guide is now ready for use.

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Model 25234

To remove the snare wire, push the finger grip forward until the loop is visible. Unscrew and remove the HF adaptor and pull the loop out of the guide tube toward proximal.

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
To insert the loop, push the finger grip forward as far as it will go and insert the new loop from distal into the guide tube as far as it will go and locate it in its seat on the grip. Screw the HF adaptor into the socket ②.

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
Machine cleaning

Make certain all parts are positioned securely to prevent damage during cleaning. The sheath tube must be rinsed out on the inside. The LUER lock irrigation port also permits easy attachment to a cleaning machine.

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Adult and pediatric proctoscopes


CAUTION: When the obturator is removed from and inserted into the proctoscope tube, the illumination rod must only be engaged in the first stage. If it is engaged in the second stage, the rod will project into the lumen and may be damaged by the obturator.

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
Insert tube 24970:

The insert tube - which takes the place of an obturator - is inserted into the proctoscope tube with the U-shape guide pointing downward.

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The insert tube is fixed by turning the set screw counterclockwise. To release the insert tube, rotate the set screw clockwise.

NOTE: An illumination rod can also be fixed in the second stage when an insert tube is inserted.

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Cleaning

CAUTION: The insert tube must not be cleaned in an ultrasound bath.

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Machine cleaning

Make certain all parts are positioned securely to prevent damage during cleaning. The sheath tube must be rinsed out on the inside.

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Anal retractor and speculum

To change the valves, pull the lock toward proximal and remove the valves toward distal.

Insert new valves into the holder and press in until they audibly snap into place.

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Suction ligature instrument

Assembly


WARNING: The knurled nuts for fastening the plunger must be between the fingers when assembling the handle, to allow the plunger to be guided. This also avoids the danger of an injury to the fingers on the handle sleeve.

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
Carefully screw in the knurled nuts up to the limit.

Place the ligature attachment onto the grip so that the 0 markings are opposite each other, and secure with the locking lever.

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
Dismantle the suction handle before cleaning and sterilizing: screw off the knurled ring and pull off the plunger.

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Remove tissue remnants with a cleaning brush. Rinse out the sheath tubes through the irrigation part. Blow dry with sterile compressed air.

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

Video operating rectoscope

For cleaning, remove the cap or attachment.

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
Remove and inspect the silicone lip valves and, if necessary, replace with new ones.

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Remove tissue remnants with a cleaning brush. Brush out and rinse out sheath tubes and parts. Blow dry with sterile compressed air.

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Inspect O-ring and, if necessary, replace with a new one.

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Operating/storage conditions

	Temperature	Rel. humidity
Storage	-20°C ... +60°C	10% ... 90%
Operation	+10°C ... +40°C	30% ... 70%

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Disposal

Standard batteries and rechargeable battery inserts may contain toxic substances, and must be disposed of in an environmentally friendly way.

To dispose of the instrument set, no special measures are necessary.

National laws and regulations must be observed.

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Directive compliance

The medical products mentioned bear the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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Limitation of liability

As the supplier of these medical products we only regard ourselves responsible for the safety, reliability and performance of the products if assembly, system expansions, readjustments, modifications, or repairs have been performed by persons authorized by KARL STORZ, and the product has been used in accordance with its operating instructions at all times.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts.

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Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective instrument.

For rigid replacement telescopes you receive a warranty of 3 years, for fiberoscopes and instruments 2 years

For videoscopes and units, individual repair is necessary.

To bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

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WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:

Company stamp/signature:

KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postfach 230, 78503 Tuttlingen, Germany
Phone: +49 (0)7461 708-0, Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de

KARL STORZ Endoskope Berlin GmbH
Ohlauer Straße 43
10999 Berlin, Germany
Phone: +49 (0)30 30 69090, Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latina-America, Inc.
815 N W 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopy
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326 , Col Chapultepec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 5250-5607, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
R. Joaquim Floriano, 413, - 20º andar -
Itaim Bibi, São Paulo - SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufralegui 627 6° Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olservei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsättrågränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4 kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytten 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopleer.be

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopiija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Gounari 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd *
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

ООО КАРЛ СТОРЗ Эндоскопия – ВОСТОК
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ Endoskope
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Sary-Arka str., 6, BC "Arman", off. 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

Karl Storz Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd,
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 811
17-1 Jeonga-dong
Seongnam, 463-847, Korea
Phone: +82-70-3450-7474, Fax: +82-70-8277-3299

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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**Important information
for users of KARL STORZ
equipment and instruments**

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Thank you for choosing a product by KARL STORZ, a name you can trust. This product, like all the others we produce, is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning and sterilization of the instrument. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

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KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology may occur for this reason. Therefore, no claims may be deduced from the information, figures and descriptions in this manual.

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Frame extension

21 Frame extension (49100 W)

- 22 Stop spring
- 23 Hole
- 24 Rotational axis
- 25 Locking screw
- 26 Holding opening
- 27 Set screw
- 28 Supporting surface
- 29 Hook
- 30 Notch
- 31 Stabilizer stand
- 32 Star knob
- 33 Locking device
- 34 Locking handle
- 35 Holding arms
- 36 Axis
- 37 Cross bar
- 38 Clamping lever
- 39 Set screw
- 40 Light cable connector
- 41 Light cable adaptor for light cable by Richard Wolf
- 42 Light cable adaptor for light cable by KARL STORZ
- 43 Clamping jaws
- 44 Unlock button (cf. p. 32)

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1 SCHOELLHORN multi-retractor System description

The SCHOELLHORN multi-retractor system is based on a concept whereby it is possible to perform all steps using a single, all-in-one retractor system.

All the components necessary for the retractor systems are used based on a modular principle. They can be customized to meet the individual requirements of the operation at hand.


Thus this system can, in addition to its function as a sternum retractor, also be used for visualizing the left and right internal mammal arteries, as for mitral valve exposition. Moreover it is employed for local stabilization of the heart in coronary bypass surgery during the 'off pump' procedure, as well as in the OPCAB and MIDCAB techniques. The system was designed to allow for optimal adaptation to the individual patient's size. A light source which can be inserted in-situ ensures excellent visibility.

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The SCHOELLHORN multi-retractor system consists of the following modules:


Standard retractor 49100 B

The standard retractor is the system's basic element and, thanks to its symmetrical shape, can be employed in any way you wish. The blades ⑤ are available in three different versions, each with a different height (4 cm as standard, 3 cm for thin patients and 5 cm for obese patients). The blades ⑤ offer notch-free adjustment and may thus be adapted to any site size. Therefore, the use of only one pair of blades is an option, making the retractor ideal for partial cranial or caudal sternotomy.

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IMA extension 49100 D


For the dissection of the internal mammary artery, an IMA elevator plate ⑥ was developed in conjunction with an opposing IMA support plate ⑨.

What's more, the retractor was designed to enable the surgeon to dissect the right internal mammary artery just as easily as the left internal mammary artery, thanks to a mirror-image system set-up.

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Mitral valve extension 49100 K


The mitral valve expansion is another of the retractor system's combination options for mitral valve exposition.

For this purpose, we developed mitral valve retractors ⑪ with a supporting surface length of 4 and 5 cm. The objective of this design is to achieve exposition of the mitral valve via the left atrial access, not by attaching the mitral valve retractors ⑪ into the atrial wall, but by tilting the mitral valve retractors ⑪ sideways. This improves the straight view, including of the anterior mitral ring, and the mitral valve level can be tilted as required. The mitral valve retractors ⑪ can also be bent with the retractor ratchet lever ④ and adapted to the anatomical conditions.

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OPCAB extension 49100 SD/49100 H

The OPCAB extension is used in coronary artery revascularization without a heart-lung machine. It consists of an OPCAB stabilizer arm ⑫ and a heart holder ⑬. The OPCAB stabilizer ⑫ can be fixed in any position. Its stabilizer stand ⑭ is flexible and can be adapted to the heart surface. To obtain a better grip on the epicardium, we recommend binding it with fabric or fixation bandage (Elastomull).

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
The heart holder ⑬ is used to hold the heart in position after it has been luxated with a cloth which is fixated behind the left atrium, to permit the posterior wall to be stabilized accordingly. The arms ⑭ of the heart holder ⑬ are also flexible and can be adjusted to the heart's dimensions. A fabric cover on the holding arms ⑮ improves the grip function.

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MIDCAB extension 49100 M

A further retractor system is necessary to perform a coronary bypass operation with lateral access using the internal mammary artery. This procedure is called 'MIDCAB'. It involves a MIDCAB elevator frame clamp ⑭, the MIDCAB elevator frame ⑮ with rotating pin ⑯, the MIDCAB elevator plate ⑰ and the MIDCAB opposite blade ⑱. After dissection of the internal mammary artery, the thoracotomy is held apart using a blade pair, and the coronary artery to be anastomosed is secured locally with the OPCAB stabilizer arm ⑫.

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
Fiberoptic light rod extension 49100 LS

Adequate illumination of the surgical site substantially improves both the dissection of the IMA via anterolateral access in a MIDCAB operation, as well as operating conditions for mitral valve surgery. This is provided in the form of a special flexible fiberoptic light rod ⑬, which is connected to a light source via a light cable, and attached to the retractor frame by a universal suspension holder ⑭.

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Frame extension 49100 W

If the surgeon chooses to place the OPCAB stabilizer arm ⑫ into the site caudally, e.g. in an OPCAB operation, since this does not interfere with the path of the internal mammary artery to the planned coronary anastomosis, he may fix the OPCAB stabilizer arm ⑫ caudally via the frame extension ⑪.

This frame extension ⑪ can easily be placed on the retractor frame and secured with two screws. It then provides a solid fourth side by means of which more components may be introduced to the surgical site.

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2 Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully read to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

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Safety instructions

Warnings and cautions

! WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no accidental rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

! WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

! WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

! WARNING: If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.

! WARNING: Mechanical damage may cause the flexible collar on the OPCAB stabilizer arm/the OPCAB extension and fiberoptic light rod to leak and liquids to enter. Therefore always inspect the collar for mechanical damage prior to using it.

i NOTE: Follow the application and storage specifications stated in the instruction manual.

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3 Intended use

The SCHOELLHORN multi-retractor is used to dilate the surgical site after sternotomy and thoracotomy.

Modular extensions for the basic retractor enable IMA dissection from both sides, surgical procedures on the mitral valve and operation of the coronaries using the 'off pump' technique for OPCAB and MIDCAB procedures, and also significantly improve visibility at the surgical site thanks to a light source.

Intended purpose of the 49100 H heart holder: Heart holder to complete the OPCAB extension comprising heart holder and OPCAB stabilizer arm, for use in cardiovascular surgery.

Intended purpose of the 49100 K mitral valve extension: Mitral valve extension comprising mitral valve retractor and holding device, to achieve exposition of the mitral valve, for use with the standard retractor in cardiovascular surgery.

Intended purpose of the 49100 SD OPCAB stabilizer arm: OPCAB stabilizer arm to complete the OPCAB extension comprising heart holder and OPCAB stabilizer arm, for use in cardiovascular surgery.

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- Mitralklappenchirurgie
- Bypasschirurgie
- MIDCAB-Eingriffe (Minimally Invasive Direct Coronary Artery Bypass)
- OPCAB-Eingriffe (Off Pump Coronary Artery Bypass)

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Indication: The following indications can be named as examples within the scope of the intended purpose:

- Mitral valve surgery
- Bypass surgery
- MIDCAB interventions (Minimally Invasive Direct Coronary Artery Bypass)
- OPCAB interventions (Off Pump Coronary Artery Bypass)

Contraindications: No contraindications relating directly to the product are known. The multi-retractor must not be used if, in the opinion of an experienced physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the method as such is contraindicated.

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4 User qualification

The multi-retractor may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the device/instrument.

The information given in these instructions only serves to instruct in the correct handling, cleaning and sterilization of the multi-retractor. It is not intended as an introduction to the surgical technique.

5 Safety precautions when operating the multi-retractor

The multi-retractor must be used according to the medical rules and procedures recognized for this method.

6 Safety precautions when handling the OPCAB stabilizer arm and fiberoptic light rod

Inspect the flexible collar on the OPCAB stabilizer arm ⑫ and the fiberoptic light rod ⑬ for cracks and scratches. All movable parts must move easily and show no sign of damage. Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced.

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Scope of supply

7 Scope of supply

SCHOELLHORN multi-retractor
(Set Art. No. 49100) consisting of:

Standard retractor (Set Art. No. 49100 B)

- 1x Fixed retractor frame arm (49100 BR)
- 1x Movable retractor frame arm (49100 BR)
- 1x Retractor pinion (49100 BR)
- 2x Retractor ratchet lever (49100 BS)
(2nd retractor ratchet lever is replacement)
- 4x 3 blade sizes
(height 3 cm: 49100 VK)
(height 4 cm: 49100 VM)
(height 5 cm: 49100 VL)

IMA extension (Set Art. No. 49100 D)

- 1x IMA elevator plate (49100 DE)
- 1x IMA assembly frame (49100 DH)
- 1x IMA notch lever (49100 DG)
- 1x IMA support plate (49100 DA)

Mitral valve extension (Set Art. No. 49100 K)

- 1x Holding device for mitral valve retractors (49100 KT)
- 2x 2 mitral valve retractors
(hook length 40 mm: 49100 KD and hook length 50 mm: 9100 KE)

OPCAB extension

- 1x OPCAB stabilizer arm (49100 SD)
- 1x OPCAB heart holder (49100 H)

MIDCAB extension (Set Art. No. 49100 M)

- 1x MIDCAB elevator frame clamp (49100 MG)
- 1x MIDCAB elevator frame (49100 MR)
- 1x MIDCAB elevator plate (49100 MP)
- 1x MIDCAB threaded rod (49100 MR)
- 2x MIDCAB opposite blade
(height 30 mm: 49100 MM and
height 40 mm: 49100 ML)

Fiberoptic light extension (Set Art. No. 49100 LS)

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8 Unpacking the equipment

Carefully unpack and remove the instruments and the accessory from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these may be useful if the instrument has to be transported.

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9 Initial use

WARNING: The instrument is not sterile when delivered, and must be cleaned, disinfected, and/or sterilized prior to initial use and prior to every subsequent use as well as immediately after every use.

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
10 Testing

All components of the multi-retractor must be checked for signs of damage prior to and following each surgical intervention. All components must be tested for correct functioning prior to use. All movable parts must move easily and show no sign of damage. All parts of the instrument must be inspected for correct orientation!

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11 Handling

11.1 Standard retractor

The free-moving, notch-free adjustable blades ⑤, can be slid onto the retractor frame. Their position can be locked by using the retractor in combination with the contact surface of the blade ⑤ with the tissue to be retracted. Because the blade elements can be positioned and adjusted infinitely, it is possible to distract the sternum mostly in the caudal half, thereby reducing the distraction of the cranial half and exposure of the brachial plexus. For optimum treatment of both thin and obese patients, blades ⑤ in heights 3 cm/5 cm can be used instead of those of the standard height of 4 cm.

The blades ⑤ are all 4 cm wide.

① **NOTE:** It is possible to use only one blade pair for especially small ports. This makes the retractor ideally suited for partial cranial or caudal sternotomy.

When positioning the basic retractor it is best to place the frame as cranially as possible so that the frame arms do not unnecessarily jut out toward caudal. The blade pairs should be positioned as far as possible in the middle and lower third of the sternum to protect the brachial plexus.

Thanks to the low profile blade sections ⑤, the retractor becomes inconspicuous and access to the surgical site is optimized and unimpeded. Flexible positioning enables optimal adaptation to the site and thus atraumatic distraction of the sternum.

Similarly, use of the standard retractor makes thoracotomy possible in a manner that is both flexible and gentle on tissues.

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
Assembly


① **NOTE:** The retractor pinion ③ is symmetrical and can be used from both sides.

Place the retractor pinion ③ in the hole ②④ on the movable retractor frame arm ②.

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
①

NOTE: Position the fixed retractor frame arm ① and the movable retractor frame arm ② so that the tab on the stop spring ②④ of the fixed retractor frame arm ① and on the movable retractor frame arm ② point in the same direction.

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
The movable retractor frame arm ② is now slid onto the fixed retractor frame arm ① with the attached retractor pinion ③ and then screwed on further by simultaneously turning the retractor pinion ③ until the teeth of the fixed retractor frame arm ① interlock.

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


Place the retractor ratchet lever ④ on the retractor pinion ③ and set the movable retractor frame arm ② to the required position.

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! **CAUTION:** Before mounting the blades ⑤, it is essential to make sure that the blades used ⑤ are of the same height.

The blades ⑤ can now be slid onto the retractor frames. The stop springs ②② are designed to prevent mounted blades ⑤ from slipping off the instrument.

! **CAUTION:** When positioning the blades ⑤ on the retractor, the rotational axes ②④ of the blades ⑤ must be inside and the rotatable retracting surface of the blades ⑤ facing inwards.

i **NOTE:** The blades ⑤ have a rhombus pattern to ensure a secure grip on the retracting surface in the surgical site.

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
Thus assembled, the retractor is ready to be used on the surgical site.

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
11.2 IMA extension

The design of the suspension mechanism for the IMA elevator plate (6) in conjunction with the opposite IMA support plate (9) means that the surgical site is well exposed. The IMA notch lever is used for further elevation (8). The contact surface of the IMA elevator plate (6), which has also been provided with a rhombus pattern, ensures a secure grip of the sternum.

Since the retractor ratchet lever (4) can be attached to either side of the retractor pinion (3), the IMA extension can be used inversely for dissection of both the left and right internal mammary artery.


To place the IMA extension on the surgical site, first insert and connect the IMA elevator plate (6). Then do the same with the retractor frame which has been preassembled with the IMA support plate (9) and two blades (5) attached to the fixed retractor frame arm (1) and the IMA assembly frame (7) on the movable retractor frame arm (2).

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
Assembly

To install the IMA assembly frame (7), remove the two blades (5) from the movable retractor frame arm (2).

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
CAUTION: When doing this, the suspension pin must point toward the surgical site.

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
To mount the IMA support plate (9), remove one blade (5) from the fixed retractor frame arm (1).

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Slide the IMA support plate (9) onto the fixed retractor frame arm (1).


WARNING: When doing this, the large surface must face outwards with the curvature facing upwards.

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Then slide the blade (5) back onto the retractor frame arm.

CAUTION: When positioning the blades (5) on the retractor, the rotational axes of the blades (5) must be inside and the rotatable retracting surface of the blades (5) facing inwards.

NOTE: The blades (5) have a rhombus pattern to ensure a secure grip on the retracting surface on the surgical site.

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

Insert the IMA elevator plate (6) in the surgical site and suspend it using the suspension pins in the IMA assembly frame (7).

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
Thus positioned and connected in the surgical site, the chest wall is well raised when the retractor is dilated.

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
In case where the elevation of the thoracic wall achieved by the geometry of the retractor is not yet not sufficient for preparation of the IMA, the IMA notch lever ⑧ can be used to improve the elevation further. The IMA notch lever ⑧ can be inserted into the cross-hole of the IMA assembly frame ⑦ from both sides, as required.

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
The IMA notch lever ⑧ should be introduced far enough to rest strainlessly in the slot of the IMA elevator plate ⑥.

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

... the elevation of the thoracic wall can be further increased by carefully pushing the IMA notch lever (8) upwards against the IMA elevator plate (6).

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
11.3 Mitral valve extension

The supporting surfaces (28) of the mitral valve retractors (11) can be bent using the retractor ratchet lever (4). Mitral valve retractors (11) are available in lengths of 4 and 5 cm. The mitral valve retractors (11) can be suspended from the holding device (10) according to the required position. They are freely height-adjustable. The valve exposition is facilitated by infinitely adjustable tilting of the mitral valve retractors (11). If necessary, the mitral valve retractors (11) can be inserted turned. In cases where the atrium is small, it is also possible to work with just one mitral valve retractor (11). When set up this way, the site affords visualization of the entire mitral valve area.

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
Assembly

Place the holding device (10) between two blades (5) on the side of the retractor frame to the left of the patient. To do this, ...

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
... the locking screw (25) must point away from the surgical site, while the holding openings (26) must point toward the surgical site to suspend the mitral valve retractors (11).

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
Secure the holding device (10) with the locking screw (25).

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
The set screw (27) can be used to tilt the holding device (10) to any position and thereby set the inclination angle of the mitral valve retractors (11). Before the mitral valve retractors (11) are inserted, the holding device (10) should be adjusted as far up as possible. In this way, after they have been positioned in the left atrium, the mitral valve retractors (11) can be positioned and then tensioned by tilting the holding device (10).

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The supporting surface (28) of the mitral valve retractors (11) can be bent using the retractor ratchet lever (4) and should be positioned to enable a direct view of the entire mitral valve.


! WARNING: The supporting surface (28) of the mitral valve retractors (11) should be inspected for cracks after each bending. If cracks are detected, the mitral valve retractors (11) must not be used.

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
... around 30 – 45° of the supporting surface (28) to the axis of the mitral valve retractors (11). When doing this, the supporting surface (28) should be as straight as possible and the deflection proximal. In this way, the view of the mitral valve is not obstructed by a crooked mitral valve retractor (11). If, as is most common, 2 mitral valve retractors (11) are used, it is important that both mitral valve retractors (11) are bent in the same angle, as this makes determining the mitral valve site easier.

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


By turning the hook (29) in different ways, the mitral valve retractor (11) can be repositioned on its axis in three notches (30) to positions A, B, or C. This allows the supporting surfaces (28) to be fixed in three directions.

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


After the hook (29) has been inserted into the holding device (10), the height can be set by pulling the mitral valve retractor (11). The hook (29) will keep the mitral valve retractor (11) tensioned by tilting (similarly to the drawer sign). By lifting the hook (29) (offsetting the tilt), the 'locking mechanism' is released again.

By introducing the hook (29) sideways into the holding device (10), the mitral valve retractor (11) is prevented from tilting away and the mitral valve retractors (11) will not slide across one another at atrium level.

In the case of a large atrium, it is therefore recommended to introduce the cranial mitral valve retractor (11) in the middle position (position B) and the caudal retractor into the site turned 30° toward caudal (position C). With a medium-sized atrium, both mitral valve retractors (11) should be in a straight position (position B) and with a small atrium it might be advisable to insert just one mitral valve retractor (11).

For better visualization of the mitral valve, attach the fiberoptic light rod extension as described from page 31 onward.


If the fiberoptic light rod extension is used, the mitral valve is visualized exceptionally well!

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
11.4 OPCAB extension

The local securing of the coronary artery required for off-pump surgery is carried out with the OPCAB stabilizer arm 12. In the case of revascularization of the posterior lateral wall, the heart holder 13 significantly facilitates stabilizing with the OPCAB stabilizer arm 12.

The contact of the stabilizer stand 31 and the holding arms 35 on the heart holder 13 with the epicardial surface can be improved by using a fabric cover, e.g. 'Elastomull' before sterilization or a sterile plaster during the operation.

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Assembly

OPCAB stabilizer arm

If required, attach the frame extension 21 as described on page 35.

Place the OPCAB stabilizer arm 12 on the required site of the retractor system or the frame extension 21 and secure using the locking screw 25.

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

Bend the stabilizer stands (11) according to the anatomy of the heart surface.


WARNING: The stabilizer stands (11) are flexible. After each bending of the stabilizer stands (11), they should be inspected for cracks. If cracks are detected, the OPCAB stabilizer arm (12) must not be used.

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Position the flexible collar of the OPCAB stabilizer arm (12) and the stabilizer stands (11) as required and turn the star knob (32) until it seizes.

CAUTION: To prevent damage to the OPCAB stabilizer arm (12), the star knob (32) must not be turned too hard. The star knob (32) should only be turned hard enough to secure the stabilizer stand (11) adequately.

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
Heart holder

Release the locking device (33) on the heart holder (13) and open the brackets (35).

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Bend the holding arms (35) according to the anatomy of the heart.


WARNING: The holding arms (35) are flexible. After each bending of the holding arms (35), they should be inspected for cracks. If cracks are detected the heart holder (13) must not be used.

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

Place the heart holder (13) on the required site of the retractor system and secure using the locking screw (25).

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
Lock the brackets (35) in the required position by closing the locking device (34).

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

Using the joints, position the holding arms 34 as required and screw the locking handle 35 tight.

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
The posterior lateral wall of the heart can be secured better by using the heart holder with the OPCAB stabilizer arm 12 in this way.

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11.5 MIDCAB extension


The MIDCAB extension was developed for lateral access for a coronary bypass operation using the internal mammary artery. The MIDCAB extension can be slid onto the standard retractor system.

The same OPCAB stabilizer arm ⑫ as that described for the OPCAB extension can be used for subsequent coronary anastomosis.

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
Assembly

Slide the MIDCAB elevator frame clamp ⑭ onto the movable retractor frame arm ② in such a way that the MIDCAB components can be suspended on the cranial side of the surgical site.

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
Position the MIDCAB opposite blade ⑮ on the fixed retractor frame arm ① so that the supporting surface of the MIDCAB opposite blade ⑮ faces the surgical site. Slide the MIDCAB opposite blade ⑮ onto the fixed retractor frame arm ①.

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
Introduce the MIDCAB elevator plate ⑮ from underneath into the MIDCAB elevator frame ⑮, ...

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
...and hang it on the axis (36) of the MIDCAB elevator frame (15) ...

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
... and using a little force, click it onto the axis (36).

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

Screw the MIDCAB threaded rod (17) into the MIDCAB elevator frame (15).

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
Hang the MIDCAB components on the cross bar (17) of the MIDCAB elevator frame (15) in the MIDCAB elevator frame clamp (14).

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
The MIDCAB threaded rod (17) allows for the height adjustment of the MIDCAB elevator plate (16).

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
If required, attach the fiberoptic light rod extension as described from page 31 onward.
Using the fiberoptic light rod (19) in addition improves the possibilities for dissection of the internal mammary artery (IMA).

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

11.6 Fiberoptic light rod extension

The flexible fiberoptic light rod (19), which is connected to a corresponding light source, is used for suitable illumination of the surgical site - particularly for the dissection of the IMA via an anterolateral access (MIDCAB) and in mitral valve surgery. The versatile fiberoptic light rod (19) optimizes the illumination of the surgical site in many cases.

Assembly

The fiberoptic light rod extension consists of a flexible light rod (19) with distal light outlet at a 90° angle, which can be rotated on the axis of the fiberoptic light rod (19) and a proximal connection for the light source. Thanks to the deflection of the light outlet, the surgical site may be illuminated while saving space and avoiding limiting surgical access. The rotatability of the light outlet means that the light beam can be positioned very accurately on the aligned light rod.

If required, screw the adaptor (41, 42) which matches the light cable used (see page IV) onto the light cable connector (40).

The fiberoptic light rod extension can be placed anywhere on the retractor system using the universal suspension holder (20).

NOTE: The rotatable end of the fiberoptic light rod (19) where the light outlet is located heats up, but should, when used correctly, not lead to burns even if placed in direct contact with tissue at the surgical site.

WARNING: Effect of heat! The end of the flexible light cable from the light source to the connector for the fiberoptic light rod (19) must not have any contact with tissue, as a large amount of heat is dissipated via the light cable in this area.

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
Slide the posterior clamping jaws (43) onto the pin of the universal suspension holder (29) and hold it in place.


NOTE: The anterior clamping jaw (41) should be held at the clamping lever (38) and not on the housing, so that the clamping lever (38) does not get pushed out of the housing when the (43) clamping jaws are attached opposite each other.

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Place the clamping lever (38) in an open position.

CAUTION: Pressing the unlock button (44) allows the clamping jaws (43) of the universal suspension holder (29) to be unlocked and disassembled.

NOTE: The stand of the fiberoptic light rod clamp (29) is designed in such a way that the fiberoptic light rod extension can be used in conjunction with other retractors too.


CAUTION: The universal suspension holder (29) is intended for square profiles from 5 x 15 mm to 8 x 35 mm. Regardless of this, functionality must always be tested before each use.

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
Place the universal suspension holder (20) on the retractor system and secure it by turning the locking screw (23).

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


Connect the fiberoptic light rod (19) to the light source at the proximal end using a suitable light cable and bend the fiberoptic light rod (19) into the required position after arranging the light beam.

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
Place the fiberoptic light rod (19) in the clamping jaws of the universal suspension holder (20). The fiberoptic light rod (19) with the clamping jaws can now be freely rotated on the mount. Once in the required position, close the clamping lever (21). This secures both the fiberoptic light rod (19) in the clamping jaws and the jaws on the mount.

CAUTION: Proper installation and function of the fiberoptic light rod clamp (20) must be checked before each use.

Subsequently align the light beam exactly so that the operating area is optimally illuminated.

NOTE: The rotatable end of the fiberoptic light rod where the light outlet is located heats up, but should, when used correctly, not lead to burns even if placed in direct contact with tissue at the surgical site.

WARNING: Pay attention to the instruction manual for the light source.

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

11.7 Frame extension

If the space on the opposing sides of the retractor frame is not sufficient to attach the OPCAB stabilizer arm (12), the fiberoptic light rod (13) or other components, the frame extension (21) may be used.

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
Assembly

Release the set screws (23) and the locking screws (24).

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
Place the frame extension (21) on the ends of the fixed retractor frame arm (1) and the movable retractor frame arm (2) and secure with the locking screws (24).

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Set the movable retractor frame arm (2) to the required position by turning the retractor ratchet lever (4).

CAUTION: The movable retractor frame arm (2) may only be adjusted when the set screws are loose (23).

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

Position the elements of the frame extension ②① and secure them with the two set screws ③⑨.

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A retractor frame extended in this way can be used for attaching various system components, even from the patient's caudal side.

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12 Disassembly

To dismantle the instrument, proceed in the same way described for assembly, but in the reverse order.

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13 Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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CAUTION: When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

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- Brushes: 27652, 27648 A, 27650 C
- Adaptor
- Cleaning gun: 27660

13.2 Disassembly

The medical device must be separated into its individual components as far as possible prior to cleaning and disinfection.

13.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

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



13.4.1 Manual precleaning SCHÖELLHORN standard retractor

Brushing the surfaces

Any visible contamination or heavy soiling must be removed from the surfaces of the fixed retractor frame arm, movable retractor frame arm, retractor pinion, retractor ratchet lever and blades by cleaning them under cold running water with the aid of a brush or sponge.

Brushing the lumina

In order to remove all visible contamination, the hole in the movable retractor frame arm, the holes in the blades and the slit in the retractor ratchet lever must be precleaned under cold running water with the aid of a suitable brush.

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
13.4.2 Manual precleaning IMA extension

Brushing the surfaces

Any visible contamination or heavy soiling must be removed from the surfaces of the IMA elevator plate, IMA assembly frame, IMA notch lever and IMA support plate by cleaning them under cold running water with the aid of a brush or sponge.

Brushing the lumina

In order to remove all visible contamination, the slots in the IMA assembly frame and IMA support plate must be precleaned under cold running water with the aid of a suitable brush.

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13.4.4 Manual precleaning OPCAB stabilizer arm

CAUTION: The flexible plastic sheath for the OPCAB stabilizer arm must not be damaged. Therefore cleaning must only be carried out using a soft disposable cloth or disinfected sponge.

CAUTION: The OPCAB stabilizer arm must not be cleaned using ultrasound. Dirt may gather underneath the flexible plastic sheath.

Brushing the surfaces

Any visible contamination or heavy soiling of the medical device surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

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For technical reasons, the medical device is not suitable for ultrasound treatment.

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
13.4.5 Manual precleaning OPCAB heart holder

Brushing the surfaces

Any visible contamination or heavy soiling of the OPCAB heart holder surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

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The hollow spaces of the OPCAB heart holder must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose.

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13.4.6 Manual precleaning MIDCAB extension

Brushing the surfaces

Any visible contamination or heavy soiling must be removed from the surfaces of the MIDCAB elevator frame clamp, MIDCAB elevator frame, MIDCAB elevator plate, MIDCAB threaded rod and MIDCAB opposite blade by cleaning them under cold running water with the aid of a brush or sponge.

Brushing the lumina

In order to remove all visible contamination, the slots in the MIDCAB elevator plate must be precleaned under cold running water with the aid of a suitable brush.

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13.4.8 Manual precleaning Universal suspension holder Brushing the surfaces

Any visible contamination or heavy soiling must be removed from the surfaces of the universal suspension holder and disassembled clamping jaws by cleaning them under cold running water with the aid of a brush or sponge.

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The hollow spaces of the universal suspension holder and cleaning holes must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose.

NOTE: There are two cleaning holes to be rinsed out (one from the front and one from the right).

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13.4.10 Manual precleaning Extension frame

Brushing the surfaces

Any visible contamination or heavy soiling of the extension frame surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

Auxiliary pretreatment with ultrasound

Ultrasound treatment at approx. 35 kHz and with a maximum exposure time of ten minutes can be used to aid the preliminary cleaning of very heavily soiled and/or complex medical devices. The device must then be rinsed with cold water to ensure neutralization.

13.5 Manual cleaning

The medical device must be completely immersed in a cleaning solution. Through disassembling and opening the instrument, and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and that no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with the aid of a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

13.6 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and that no air bubbles are present. Following the exposure time, the medical device must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all the surfaces, joints, openings, channels and lumina are dried completely with sterile compressed air. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose.

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13.7 Machine cleaning and disinfection

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Thermal disinfection is preferable. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

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13.8 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Glide surfaces must be lubricated locally with instrument oil.
- Dismantled medical devices must be assembled (see Chapter 11).
- Afterwards, a functional check must be carried out.

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NOTE: The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

13.9 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953).

13.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

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The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134 °C – 137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes. These procedures are only suitable for thermostable instruments.

NOTE: Greased components must be disassembled for sterilization, so that steam penetration is ensured.

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– 39360 AK

13.11 Accessories

Accessories recommended for sterilization and storage:

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14 Servicing and repair

Defective items of equipment should be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

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15 Limitation of liability

KARL STORZ GmbH & Co. KG shall be liable for safe operation, operational reliability, and performance of instruments only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that instruments have been used in accordance with its operating instructions at all times.

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16 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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17 Directive compliance

These medical products bear the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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18 Disposal

No special measures are necessary to dispose of these medical products. National laws and regulations must be observed.

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19 Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which usually ensures the prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective instrument. You receive a warranty of 1 year for telescopes and 6 months for other instruments.

For miniature and neuro-endoscopes, individual repair is necessary. To bridge the repair period, you will receive an instrument on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany, you can send instruments to be repaired directly to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries, please contact your local KARL STORZ branch or authorized dealer.

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20 Important information

To avoid the spread of infectious illnesses among hospital staff or KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 9993390, Fax: +961 1 9993391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kazakhstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O. Box 118069
Dubai – United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/719) Dang Van Ngu,
F10-Q, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renmin Road South
Chengdu, Sichuan, 610016
People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjae-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisan Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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14. Cleaning, disinfection, sterilization31	
14.1 Reprocessing Overview 1, German.....34	
14.2 Reprocessing Overview 1, English35	
14.3 Reprocessing Overview 1, Spanish.....36	
14.4 Reprocessing Overview 2, German.....37	
14.5 Reprocessing Overview 2, English38	
14.6 Reprocessing Overview 2, Spanish.....39	
14.7 Reprocessing Overview 3, German.....40	
14.8 Reprocessing Overview 3, English41	
14.9 Reprocessing Overview 3, Spanish.....42	
15. Preliminary cleaning43	
15.1 Disassembling suction valve insert45	
16. Manual reprocessing.....46	
16.1 Manual cleaning.....46	
16.2 Manual disinfection48	
17. Machine reprocessing50	
18. Inspection after reprocessing58	
18.1 Assembling suction valve insert59	
19. Gas sterilization/ plasma sterilization60	
20. Chemical sterilization.....63	
21. Storage64	
22. Accessories.....65	
23. Overview: Cleaning brushes.....68	
24. Normatives70	
24.1 Servicing and repair70	
24.2 Limitation of liability70	
24.3 Warranty70	
24.4 Directive compliance70	
25. Repair and replacement program...71	
25.1 Important information.....71	
Subsidiaries.....72	

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Important note for the users of Broncho-Fiberscopes and Intubation Fiberscopes

It is recommended to check the suitability of the product for the intended procedure prior to use.

It has been determined that Broncho-Fiberscopes and Intubation fiberscopes (see table next page) incorporate a lumen material that is only approved for use as an external surface material and not included in the STERRAD® sterilization systems FDA-cleared claims for use in flexible endoscope lumens. This could result in lower than expected levels of hydrogen peroxide within the lumen during the exposure phase of the sterilization cycle and compromise sterilization efficacy.

Based on the risk factors discussed above, the probability of a patient infection occurring from a cleaned and disinfected bronchoscope or intubation scope processed in the STERRAD® sterilization systems is defined as "remote". Since bronchoscopes and intubation scopes are semi-critical (B) devices and need not necessarily sterilization, KARL STORZ has validated many sterilization methods.

Available validated sterilization methods include:

- STERRAD® 100S with Booster,
- STERIS® System 1,
- STERIS® System 1E,
- STERIS® AMSCO® V-PRO® maX and
- ethylene oxide (EO) gas.

The instructions „Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments“, Version 1.0.0, No. 96216003DER, describe in detail the workflows and process-relevant parameters of the individual validated procedures. The procedure must be chosen taking into account applicable national requirements and in consultation with the device manufacturers.

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The methodologies for manual and automated cleaning and disinfection are already described in details in this instruction manual. When selecting the disinfectant the range of action should be taken into account with regard to microorganisms. Use only chemicals approved by KARL STORZ for processing. For a complete list, please visit www.karlstorz.com. When preparing and using solution, strictly observe the chemical manufacturer's instructions on concentration, temperature, exposure time, and shelf life.

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Fiberscopes Series: 110xx ..., 111xx ..., 112xx ..., 113xx ... *)

*) x = number; ... = letters
e. g.: 11301 BN or 11301 BNP

NOTE: Separate instruction manuals are available for flexible videoscopes; they are therefore not described in this brochure.

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
Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of fiberscopes manufactured by KARL STORZ. The fiberscopes are sensitive precision instruments which require particular care during reprocessing. To prevent damage during reprocessing, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ fiberscopes. For this reason, read these instructions carefully and keep them as a ready reference.

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DCI®-Fiberskope (Direct-Coupled-Interface)

This instruction manual describes three fiberoptic versions, distinguished from one another by the imaging system. If a section refers exclusively to one of these versions, it is thus appropriately labeled.

Telescope

Telescope with special filters for PDD/AF

DCI® fiberoptic (Direct Coupled Interface)

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


KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope for supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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2. Controls


Connectors, accessories and their uses






- 495 F Light cable adaptor for Wolf
- 495 G Light cable adaptor for KARL STORZ, Olympus, Winter & Ibe
- 11025 E Pressure compensation cap (red)
- 11301 CA Lip valve, disposable version, pack of 20
- 13242 XA Adaptor for leakage test in cleaning and disinfection machines (green)
- 29100 Plug for instrument channel, pack of 10, reusable
- 11301 CB Suction valve insert, reusable version
- 5917800 Seal ring (suction valve insert)
- 6885391 Cap (suction valve insert)
- 5923110 Spring (suction valve insert)
- 6927791 Housing (suction valve insert)
- 8567290 Plunger (suction valve insert)
- 11301 CE Suction valve insert, disposable version, pack of 20
- 11301 D1 LED battery-powered light source with fine screw thread
- 11301 D3 LED battery-powered source with quick-release thread
- 6927691 Three-way double stopcock adaptor with LUER-Lock
- 6011590 LUER cap
- 600007 Barbed tube, thin, with LUER-Lock

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
20260030 DCI® adaptor, for connecting endoscopes with DIN eyepiece to DCI® camera heads

20260031 DCI® eyepiece adaptor for connecting DCI® telescopes to camera heads with mounts for DIN eyepieces

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- ① Focusing ring (if available)
- ② Suction valve insert (if available)
Press = suction activated
Release = suction deactivated
- ③ Valve nozzle (if available)
- ④ Up/down control
- ⑤ Valve for pressure compensation
and leakage tester connection
- ⑥ Suction outlet
- ⑦ Light inlet piece
- ⑧ Instrument channel
- ⑨ Objective lens
- ⑩ Light outlet
- ⑪ Instrument in instrument channel

- ⑫ PDD filter changer:
ALA for special tumor marker
based on aminolevulinic acid
W for white light
AF autofluorescence mode

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3. Intended Use

Fiberscopes are used in diagnostic and surgical endoscopy. They must only be used by appropriately trained and qualified technical staff.

Indications: Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in various medical fields.

Contraindications: No contraindications are known which directly relate to the product. The fibroscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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3.1 Safety precautions when using the fiberscopes

Fiberscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

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4. User qualification

KARL STORZ fiberscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique.

The information given in these instructions only serves to instruct in the correct handling, cleaning, disinfection and sterilization of the fiberscopes. It is not intended as an introduction to the endoscopic technique.

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5. Warnings and cautions

Please read this manual and follow its instructions carefully. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the text should be carefully reviewed to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

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CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the product.

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WARNING: Before using the fiberscope, read the following safety instructions carefully, to avoid putting your patients, personnel or yourself at risk.

WARNING: Inspect the fiberscope prior to each surgical procedure to ensure that it is functioning correctly and that it has been correctly cleaned, disinfected and (where applicable) sterilized.

WARNING: This equipment may not be used during the discharge of a cardiac defibrillator. Remove all endoscopic equipment from patient contact and from the surgical site.

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Warnings and cautions

WARNING: Verify the sound condition of the fibroscope as well as of accessories used in combination before every application. Damaged fiberscopes or damaged accessories must not be used.

WARNING: Combinations of medical devices are only then assured to be safe if

- they are identified as such in the respective instruction manual or
- the intended purpose and the interface specifications of the devices used in combination permit this.

Pay careful attention to the instructions and interface specifications of medical products used in combination.

WARNING: Lengthy direct contact between the endoscope tip and surrounding tissue may cause tissue damage due to the potentially high temperatures on the tip of the endoscope and must thus be avoided.

WARNING: The fibroscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

WARNING: Remove the red pressure compensation cap from the vent port

- for every endoscopic intervention
- before immersion in liquid/chemical sterilization.

WARNING: Before the fibroscope is removed from the patient, the locking mechanism (5) must be released and the distal tip of the fibroscope put in the neutral, non-angled position. Otherwise the instrument may be damaged or the patient injured.

WARNING: If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the fibroscope in the neutral, non-angled position and remove the fibroscope slowly and carefully from the patient.

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CAUTION: Do not steam sterilize (autoclave) the fiberscopes.

CAUTION: The maximum permissible temperature for cleaning, disinfection, sterilization, transport and storage is 65 °C.

CAUTION: Never immerse fiberscopes in normal saline solution, since even only brief contact can lead to pitting and corrosion.

CAUTION: Material damage: if the outer casing of the fiberscope is cleaned with alcohol, this will damage the material.

CAUTION: Fiberscopes must not be cleaned in an ultrasound bath under any circumstances.

CAUTION: Place the red pressure compensation cap on the vent port before

- gas sterilization/plasma sterilization
- shipping

CAUTION: Only use KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

CAUTION: Only use guide probes with flexible ends. Other types of guide probe may cause considerable damage to fiberscopes.

CAUTION: Make certain to read the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. Methods of cleaning, disinfection and sterilization are explained in detail there.

NOTE: Any damage to the fiberscope resulting from incorrect operation is not covered by the guarantee.

NOTE: Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

NOTE: If the instrument is in the working channel, its additional rigidity may limit angling.

NOTE: Separate instruction manuals are available for flexible videoendoscopes; they are therefore not described in this brochure.

NOTE: Federal law restricts this device to sale by or on the order of a physician (USA only).

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6. Unpacking

Check for missing items and evidence of shipping damage. File any complaints to KARL STORZ or the supplier immediately.

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6.1 Basic equipment

	Art. no.
1 Fiberscope	
1 Suction valve insert, disposable version, pack of 20 inserts	11301 CE
1 Leakage tester	13242 XL
1 Pressure compensation cap	11025 E
1 Plug for instrument channel, pack of 10	29100
2 Adaptor for cleaning the working channel in cleaning machines (green)	11301 CD
1 Cleaning brush for 2.8 x 50 and for 2.8 x 70	27651 A*)
for 3.7 x 65	27651 AL*)
for 5.0 x 65	27651 B
for 5.6 x 65	27651 C
for 6.4 x 65	27651 C
1 Biopsy forceps (only bronchofiberscopes)	11001 KL
1 Grasping forceps (only bronchofiberscopes)	11002 KS
1 Three-way double stopcock adaptor with LUER-Lock for O ₂ connection (only intubation fiberscopes)	6927691
1 LIPP tube holder (only intubation fiberscopes)	11301 CF
1 Mouth piece (only gastro-, intubation and bronchofiberscopes)	13272
*) For the fiberscopes 11278 A/AI/AU/AIU use the cleaning brush 27651 A (27651 AL is not suitable for these fiberscopes).	

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6.1.1 Accessories

- | | |
|--|----------------------|
| 1 Suction valve insert,
reusable version | 11301 CB |
| 1 Battery light source | 11301 D |
| 1 Adaptor for leakage test in
cleaning and disinfection
machines (green) | 13242 XA |
| 4 Cleaning adaptors for
cleaning machines | 11301 GG2/FF2/HH/KK2 |

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7. Initial operation

WARNING: The fiberoscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

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7.1 Fitting the suction valve

Place the suction valve insert ② on the valve nozzle ③ and turn until the projections on the side of the insert lie over the nozzle recesses.

Push gently until the suction valve insert clicks into the valve nozzle.

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7.2 Fitting the tube holder (only intubation fiberoscopes)

Push the tube holder onto the fiberoscope from distal and secure it by turning the two parts in opposite directions.

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7.3 Fitting the lip valve (only bronchofiberscopes)

Place the disposable lip valve onto the instrument channel ⑧ and tighten gently by turning it in a clockwise direction (plastic part).

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8. Handling the fiberscopes

The optical section of KARL STORZ fiberscopes consists of a fiberoptic bundle in a flexible and watertight casing, a light guide, as well as objective and eyepiece lenses. Like all optical systems, these parts are very delicate. For this reason, careful handling of the fiberscopes is essential for a long service life and a consistent quality. Therefore be aware of the following points:

- Never kink or tightly coil the fiberscopes and light guides, so that no optical fibers break off. It is therefore best to suspend the fiberscopes for storage or keep in large containers.
- Do not jerk the up/down control ④ or move forcefully.
- Do not insert the instruments or cleaning brushes forcefully through the channels. This may cause damage to the casing, especially if the tip of the instrument is deflected.
- Carefully protect the objective lens ⑨ at the distal end from damage. Do not hit. Never remove soiling with sharp-edged objects.
- Never store or disinfect the fiberscope together with sharp-edged objects to prevent damage to the exterior.
- Never place heavy objects on top of the fiberscope.

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NOTE: Every fiberscope supplied by KARL STORZ includes a photo showing the condition of the fiber bundle upon delivery. A copy of this photo is kept on file in our quality control office. The wedge-shaped marking (at 12 o'clock) serves as a guide for the user.

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9. Using the fiberscopes

The following instructions are recommended for the proper operation of the fiberscope and accessories. This is not a reference for endoscopic techniques.

Consult the appropriate medical literature for information on endoscopic techniques.

You can find accessories for your fiberscope (catheters, flexible scissors and forceps, guide wires, stone baskets, lithotripter probes, laser fibers, electrodes etc.) in the appropriate KARL STORZ catalog.

WARNING: Always adjust the light source to the minimum illumination intensity necessary to achieve optimum visualization of the operating field. The higher the intensity setting of the light source, the greater the heat energy that will be generated at the tip of the fiberscope.

WARNING: Never lay the end of the fiber optic cable or the distal tip of the fiberscope on or under a sterile cover cloth. The light intensity may lead to burns on the patient and/or sterile cover cloth.

- When connecting and operating light sources and suction/irrigation equipment, be certain to read the appropriate instruction manuals for the equipment. The irrigation equipment recommended by KARL STORZ is set to a maximum pressure of 54 kPa (400 mmHg/0,54 bar/7.8 psi).
- Endoscopic examinations must be carried out under aseptic conditions.
- When introducing the fiberscope into the patient, the distal tip of the fiberscope must be in a straight, nondeflected position (control lock-in can be felt and heard). The fiberscope must be handled according to the recognized medical rules and procedures for endoscopic methods.

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- If no instrument is inserted, close the instrument channel proximally with a plug (Art. no. 29100).
- If an instrument is to be inserted, position the disposable lip valve (Art. no. 11301 CA) to seal off the instrument channel proximally.
- Instruments especially designed for this may be inserted through the instrument channel of the fiberscopes. Make certain the accessories have a suitable length and diameter. Unsuitable accessories may damage the fiberscopes. For questions on accessories, please contact KARL STORZ directly or an authorized dealer.

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! **NOTE:** When the lip valve is first used, the instrument is subjected to increased resistance. This resistance will decrease the more often the lip valve is used.

- When inserting an instrument into the instrument channel, the fiberscope must be in a straight, nondeflected position.

! **CAUTION:** When inserting instruments, increasing resistance is felt with increasing deflection. Forced introduction of instruments can lead to damage to the fiberscope or the instrument channel.

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CAUTION: If undue resistance is encountered when inserting an instrument, do not attempt to force the instrument. Withdraw the instrument and apply a water soluble lubricant to the tip of the instrument and attempt to insert again. Prior to removal of the instrument from the patient, the distal tip must be in a straight position.

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- When using forceps, scissors or other accessories that can be opened or closed, ensure that the instrument is completely closed prior to insertion into the instrument channel. Do not attempt to open the jaws until the distal tip is completely free of the channel. Opening the jaws inside the instrument channel can lead to damage to the fiberscope.

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WARNING: When using HF electrodes or laser probes, make sure that the active parts are clear of the instrument channel and in full view of the operating surgeon before power is applied. Activation within the instrument channel may damage the fiberscope and/or injure the patient.

WARNING: Make certain to read the instructions for the laser unit and/or HF generator.

CAUTION: If accessories with sharp tips, such as EHL probes or laser probes, are inserted into the instrument channel of a deflected fiberscope, this can damage the fiberscope.

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CAUTION: When using a stone basket, the stone must not be retracted through the instrument channel but instead the entire fibroscope must be removed. Once the stone is removed, the stone basket may then be taken out of the instrument channel.

CAUTION: The fiberscopes must be cleaned and then reprocessed immediately after the intervention in order to prevent incrustations. Dried-on incrustations may jeopardize successful decontamination and disinfection.

- A leakage test must be carried out after every intervention, to prevent damage during reprocessing.

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CAUTION: Use only KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

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NOTE: If the instrument is in the working channel, its additional rigidity may limit maneuverability.

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
10. Preparation of the fiberscope

Before using the fiberscope on a patient, exterior damage and operational problems must be eliminated.

Check the outer sheath of the fiberscope. Be especially aware of visible damage to the casing, particularly in the area of the fiberscope tip.

Look through the eyepiece and make certain that you have a clear image. The grid (Moiré effect) should be easily recognizable. If the image is not clear but cloudy, clean the eyepiece and lens with a cotton applicator soaked in alcohol, maintaining hygienic conditions at all times.

DCI® fiberscope image quality can be controlled on screen. If necessary, clean image waveguides, light cables and eyepieces as described above.

If it is still cloudy, it is possible that the objective lens is scratched or that there is moisture inside the telescope. The fiberscope must not be used in this case, but instead sent in for repairs.

! **CAUTION:** Before every application on the patient, check the instrument for leaks.

Check the fiberscope for leaks before every application on a patient.

Make certain to read the specifications on page 43 of this instruction manual.

10.1 Connect these supply lines

- If necessary, remove the adaptor 495 G or 495 F from the light inlet piece and connect the fiberoptic light cable or battery-powered light source to the light inlet piece ⑦.

i **NOTE:** Not necessary for battery-powered light sources with coarse threads (11301 D3).

! **WARNING:** The instruction manuals of the supply units used must be observed.

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- suction to suction outlet ⑥
- or O₂ supply (under anesthesia) to instrument channel ⑧ via the appropriate adaptor

CAUTION: The instruction manuals of the supply units used must be observed.


Switch on the respective supply units.

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Insert the DCI® fiberscope into the DCI® camera head until it locks into place.

Insert the light cable (Art. no. 495 DV) into the DCI® camera head as far as it will go.

WARNING: The instruction manuals of the supply units used must be observed.

Connect the light cable and the video signal cable up to the appropriate supply units.

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
10.2 Controlling the distal tip

The up/down movement of the distal tip on the fiberscope can be adjusted at the control ④: move smoothly and slowly.

Some fiberscopes have:

- contrapositive deflection direction (urology);
- a locking mechanism ⑮ for securing the fiberscope tip (gynecology).

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Models with locking mechanism ⑮:

To secure the deflection position of the distal tip, turn the locking mechanism ⑮ slowly to distal until resistance can be felt.

CAUTION: Forcing the locking mechanism ⑮ further in the same direction may damage the locking mechanism ⑮.

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! WARNING: Before the fiberscope is removed from the patient, the locking mechanism (15) must be released and the distal tip of the fiberscope put in the neutral, non-deflected position. Otherwise the instrument may be damaged or the patient injured.

To release the locking mechanism (15) from the secure position, turn it slowly to proximal.

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10.3 Suction

Press the suction valve (2) to start suction (when the suction pump has been connected and switched on).

! CAUTION: There should be no instrument in the suction channel during suction. The instrument channel inlet should be closed so that no secondary air is drawn in.

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10.4 Insufflation


Under anesthesia, connect the O₂ hose to the instrument channel:

- using a barbed tube (Art. no. 600007) on the instrument channel or
- if you are using the three-way double stopcock adapter with LUER-Lock (Art. no. 6927691), attach the hose to the tube holder on the angled part of the adapter. This ensures that the straight-line passage of the adapter remains free for instrument insertion.

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10.5 Focusing

The image is focused with the focusing ring (1).

i NOTE: When connecting photo or video systems recommended by KARL STORZ, the fiberscope focusing ring must be at the zero position (mark).


DCI® fiberscopes are focused using the adjustment wheel of the DCI® camera (see page 28).

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10.6 Instrument channel


Flexible instruments can be inserted through the instrument channel ⑧ as required.

If the instrument channel is not being used: replace the sealing cap.

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10.7 PDD version

Fiberscopes intended for use in Photodynamic Diagnosis (PDD) or autofluorescence (AF) are illuminated with a special diagnosis light (D-Light C). Adjust the PDD filter changer ⑫ for special filters designed for this light and possibly for marker substances:

ALA	for a special tumor marker based on aminolevulinic acid
W	for white light
AF	autofluorescence mode

① **NOTE:** If the filter changer is set to 'ALA' or 'AF', a yellowish image is produced when viewing through the integrated special filter.

① **NOTE:** No filter changer is available on PDD fiberscope art. no. 11272 CI. The filter (ALA) is permanently installed.

For optimal illumination, use the fluid light cable 495 FP (2.5 m) or 495 FO (1.8 m) between the light source and PDD fiberscope.

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
11. DCI® fiberscope with camera

11.1 Assembly

Slide the DCI® into the snap-in DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

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Connecting the special light cable 495 DC or 495 DV to the DCI® camera head:


Insert the light cable plug into the DCI® camera head until it locks into place.

While doing this, hold the light cable by the handle piece. Never pull on the cable!

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
NOTE: The entrance and exit faces of the light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Clean the surfaces with a cotton applicator soaked in cleaning solution, then disinfect it with a cotton applicator soaked in 70 % isopropyl alcohol and dry.

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DCI® fiberscopes are focused using the adjustment wheel on the DCI® camera.



The camera control unit is controlled using the two programmable control keys (silver, blue) on the DCI® camera head.

NOTE: Further controls are described in the instruction manual of the camera control unit used.

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11.2 Disassembly


Hold the DCI® camera and the DCI® fiberscope still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released.

Pull the light cable out of the DCI® camera head.

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

12. DIN eyepiece on DCI® fiberscope

The DCI® eyepiece adaptor 20260031 ensures the connection of DCI® telescopes to camera heads with DIN eyepiece mounts.

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12.1 Assembly

Slide the DCI® fiberscope into the 'snap-in' DCI® coupling mechanism of the DCI® eyepiece adaptor until the coupling snaps in automatically.

Connect a light cable or a battery-powered light source.

WARNING: The instruction manuals of the supply units used must be observed.

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12.2 Disassembly


Remove the light cable or battery-powered light source.

Hold the DCI® eyepiece adaptor and the DCI® fiberscope. Press the release button on the DCI® eyepiece adaptor. The 'snap-in' DCI® coupling mechanism is released.

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13. DCI® camera on DIN eyepiece

13.1 Assembly

Put the focusing ring (if available) of the fiberscope in the zero position (mark).

Slide the DCI® adaptor 20260030 into the 'snap-in' DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

Connect the standard light cable at the light inlet piece ⑦ of the fiberscope.

NOTE: The special light cable 495 DC or 495 DV cannot be used and may have to be removed from the DCI® camera.

WARNING: The instruction manuals of the supply units used must be observed.

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Opening the instrument coupling: turn the outer ring against the spring force, insert the fiberscope DIN eyepiece and let go of the outer ring. The instrument coupling closes automatically as a result of the spring force.

CAUTION: No further tightening of the outer ring to position the eyepiece is required. This may, among other things, cause damage to the instrument coupler!

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13.2 Disassembly

Hold the DCI® camera and the fiberscope still. Open the instrument coupler for the DCI® adaptor and remove the fiberscope.

Hold the DCI® camera and the DCI® adaptor still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released. Remove the DCI® adaptor.

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14. Cleaning, disinfection, sterilization

Cleaning, disinfection, care and sterilization

! **CAUTION:** The manual 'Cleaning, sterilization and care of KARL STORZ instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail.

! **CAUTION:** For cleaning and disinfection, use a product from the KARL STORZ approved products list. The current list of approved products can be found on our website (www.karlstorz.com).

! **CAUTION:** When preparing and using the solutions, follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion and incorrect concentrations may result in damage. Bear in mind the microbiological range of action of the chemicals used.

! **WARNING:** The stated cleaning, disinfection and, where necessary, sterilization methods alone do not guarantee proper reprocessing or sterility. This can only be achieved if the staff follows recognized and validated reprocessing methods.

! **WARNING:** For all cleaning and disinfection of contaminated endoscopes and instruments, observe the guidelines of the Employers' Liability Insurance Association and equivalent organizations.

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Cleaning, disinfection, sterilization

WARNING: Remember that sterilization can only be performed successfully if

- all organic and chemical residue on the instruments has been eliminated. Sterilization is only effective on clean, disinfected surfaces!
- work proceeds according to recognized and validated methods.
- any deviation from recommended procedures must be validated by the user.
- the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.

CAUTION: Never clean fiberscopes with a metal brush.

CAUTION: Constantly changing between different reprocessing methods causes particular wear on the materials and must therefore be avoided. A once-only change in methods however, e.g., when acquiring new cleaning equipment, is perfectly safe.

CAUTION: Fiberscopes must not be cleaned in an ultrasonic bath.

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CAUTION: Never immerse fiberscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.

CAUTION: Final rinsing during manual reprocessing of fiberscopes should be carried out using microbiologically pure/sterile water.

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Reprocessing Overview 1										
Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine	Cleaning machine	Chemical	Steam	Gas		Plasma	Chemical	
		60 °C	93 °C		134 °C, fractionated	EO	FO	STERRAD	STERIS	
110xx ..., 111xx ..., 112xx ..., 113xx ... Fiberscopes	✓	✓	-	✓	-	✓	✓	✓	✓	
495 F/G; Light cable daptor	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11025 E; Pressure compensation cap	-	-	-	-	-	✓	✓	✓	-	
11301 CD; Adaptor for cleaning working channel in cleaning machines	✓	✓	✓	✓	-	✓	✓	✓	✓	
11301 CF; LIPP tube holder	¹ ✓	¹ ✓	¹ ✓	✓	✓	✓	✓	✓	✓	¹) dismantle for cleaning
11301 D1, 11301 D3; Battery light source (accessory)	⁵ ✓	-	-	⁵ ✓	-	✓	✓	✓	✓	⁵) Battery compartment with closed cover; liquid <15 cm deep; no pressure impact; rinse only with unpressurized liquids; do not dry with compressed air
11301 FF2; Adaptor for leakage test in ETD2/ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	

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14.5 Reprocessing Overview 2, English
Preliminary cleaning

- Procedure not appropriate!
- ✓ Procedure appropriate!

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Reprocessing Overview 2										
Product	Manual	Cleaning		Disinfection	Sterilization					Remarks
		Cleaning machine	Cleaning machine	Chemical	Steam	Gas		Plasma	Chemical	
		60 °C	93 °C		134 °C, 3 minutes	EO	FO	STERRAD	STERIS	
11301 GG2; Adaptor for cleaning irrigation and working channel in ETD2/ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11301 HH; Adaptor for leakage test in BHT cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11301 KK2; Adaptor for cleaning working channel in ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
13242 XA; Adaptor for leakage test in cleaning und disinfection machines	2✓)	2✓)	2✓)	✓	✓	✓	✓	2✓)	✓	2) only PH neutral cleaners 3) change in color possible
13242 XL; Leakage tester	-	-	-	-	-	-	-	-	-	
13272; Mouth piece	✓	✓	✓	✓	✓	✓	✓	✓	✓	
29100; Plug for instrument channel	✓	✓	✓	✓	✓	✓	✓	✓	-	
600007; Barbed tube, thin, with LUER-Lock	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6011590; LUER cap	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6927691; Three-way double stopcock adaptor with LUER-Lock	✓	✓	✓	✓	✓	✓	✓	✓	✓	Disassemble stopcocks before cleaning and sterilizing

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Reprocessing Overview 3										
Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine	Cleaning machine	Chemical	Steam	Gas		Plasma	Chemical	
		60 °C	93 °C		134 °C, fractionated	EO	FO	STERAD	STERIS	
11301 CB; Suction valve insert, reusable version	4✓	4✓	4✓	✓	✓	✓	✓	✓	✓	1 dismantle for cleaning
5917800; Seal ring (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
5923110; Spring (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6885391; Cap (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6927791; Housing (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
8567290; Plunger (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
20260030; DCI* adaptor	✓	✓	-	✓	-	✓	✓	✓	✓	
20260031; DCI* eyepiece adaptor	✓	✓	-	✓	-	✓	✓	✓	✓	

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15. Preliminary cleaning

Remove any heavy soiling (secretions, etc.) from the outer casing of the fibroscope immediately after use. Use a disposable cloth moistened with cleaning solution (no solutions predominantly containing alcohol). Submerge the distal end in a container with cleaning solution and activate the suction to clean the working channel.

CAUTION: Before immersing the instrument in cleaning and disinfectant solution, check that it is leakproof, as otherwise penetration of liquid may cause damage.

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
Connect the leakage tester (art. no. 13242 XL) to the valve for pressure compensation ⑤ and build up pressure until the pointer is in the blue marked range. Watch the pointer for approx. 30 seconds: the pressure must not drop during this time. Press the red button on the pressure valve and release the pressure. Remove the test unit from the fibroscope.

NOTE: During inflation, make certain that the red pressure relief button is completely closed.

If there is a loss in pressure, the fibroscope has a leak and must neither be used nor immersed in liquid. In this event, send the fibroscope in for repairs.

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WARNING: Risk of infection: In order to prevent infection, it is strictly forbidden to ship contaminated medical devices! All medical devices must be decontaminated on-site, to avoid contact and aerogenous infections among staff. To eliminate the risk of infection, before the fiberoptic is dispatched it should be wiped-down disinfected and sterilized in the NT procedure with the pressure compensation cap in place. We reserve the right to return contaminated instruments/equipment to the sender.

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- Remove all instruments from the instrument channel ⑧.
- Remove adapter and other removable parts from the light inlet piece ⑦ and instrument channel ⑧.
- Remove suction valve insert ② from the valve nozzle ③: press the projections on the side together and pull out the suction valve insert.
- Dispose of disposable parts removed. If necessary, disassemble other parts further and reprocess separately.

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15.1 Disassembling suction valve insert


- ① **NOTE:** Only the multi-use suction valve insert (silver-coloured cap, metal plunger) is suitable for reprocessing.

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
16. Manual reprocessing

NOTE: Before cleaning, remove both adaptors from the light inlet piece ⑦ to permit free access to the fiber surface.

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16.1 Manual cleaning

Submerge the instrument in a cleaning solution approved by KARL STORZ and clean it using a cleaning agent (brush or cleaning gun, Art. no. 27660). Avoid splashing with the cleaning agent.

CAUTION: Pass the brush through the channel in only one direction and pull back only after the brush head has completely left the distal end. The brush must not be pulled back and forth within the channel, since this can lead to damage.

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Remove the brush from the channel carefully. Then thoroughly rinse the channels repeatedly using cleaning and disinfection agents. If required, add a cleaning agent booster.

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Detach the adaptor, clean the stopcocks and grease (see the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions).

The adaptor 11301 CD can be used for manual rinsing of the suction/irrigation channel.

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
Rinse to remove the cleaning solution.

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Clean the eyepiece and objective lens with a soft cloth, sponge or cotton swab, soaked with 70 % alcohol.

CAUTION: Cleaning with the special cleaning paste should only be carried out if the image as viewed through the endoscope is visibly cloudy or if there are visible deposits on the optical surfaces, it should certainly not be used routinely during reprocessing.

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Stubborn soiling and deposits may be removed with cleaning paste (art. no. 27661). To do this, a small amount of cleaning paste is placed on a clean, moistened cleaning swab and distributed over the glass surface. Through polishing, even stubborn deposits can be removed.

NOTE: The instrument must be reprocessed again following treatment with the cleaning paste!

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Finally, dry the surfaces of the fibers and the optical end faces with a cotton-tipped applicator soaked in 70 % alcohol.

NOTE: The fiber surface in the light inlet piece ⑦ must be dried off especially thoroughly with alcohol. Residue from disinfectants and detergents in the light inlet piece could bake on if the fiber optic light cable is connected and severely impair light transmission.

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16.2 Manual disinfection

All parts of the KARL STORZ fiberscope can be disinfected. For this purpose, materials approved by KARL STORZ and recommended for use with endoscopes should be used. A list of the agents recommended and approved by KARL STORZ for manual cleaning/disinfection can be found on the Internet at www.karlstorz.com/Hygiene.

For disinfection, we recommend a plastic basin (Art. no. 27647 V), in order to avoid scratches on the optical end faces in particular. A plastic basin also prevents electrolytic corrosion, which can occur if different metals are present in a solution.

For disinfection, submerge the fiberscope completely. Ensure that the lumina of the channels and sheaths are completely filled with the solution. There must be no more air (air bubbles) in the channels. Make sure that all the stopcocks, the valve insert and the valve cap are open or have been disassembled.
Immerse the valve insert, adaptors etc. separately.

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Rinse out the instrument channel ⑧ with a syringe filled with disinfectant solution. Follow the chemical manufacturers' instructions for concentration and exposure times.

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For each rinse, use new microbiologically pure/sterile water. It is absolutely necessary to rinse the fiberscope with microbiologically pure/sterile water in order to prevent potential toxic impact from any remaining chemical residue, which could cause the instrument to become conductive. This is also required to optimize the results of sterilization and avoid damage to the surface.

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
Dry the outside of the instruments with a lint-free, cloth and/or sterile compressed air.

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Remove any water in the instrument and suction channels using sterile compressed air.

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CAUTION: The pressure must not exceed 50 kPa (0.5 bar/7.25 psi). The cleaning gun might therefore have to be operated with a pressure reducer.

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NOTE: If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberoptic scopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.

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NOTE: For flexible fiberoptic scopes without channels, a hanging storage position is not essential.

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17. Machine reprocessing

After appropriate preliminary cleaning

After the appropriate preliminary cleaning (see p. 43), all KARL STORZ fiberscopes may be connected to devices for semi-automatic or automatic reprocessing. The choice of reprocessing programs must be agreed with the machine manufacturer. There are usually special programs available for fiberscopes which are suitable for cleaning and which are gentle on materials. For machine reprocessing, make sure that the fiberscopes are not heated above 65 °C.

CAUTION: Check the fiberscope for leaks before every reprocessing in semi- or fully automatic equipment.

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Leakage tester connection on washer/disinfector:

- ETD3/ETD2 washers: Adaptor for leakage test (Art. no. 11301 FF2, see pp. 52/55)
- BHT washer/disinfector: Adaptor for leakage test (Art. no. 11301 HH)
- Belimed washer/disinfector: Adaptor for leakage test (Art. no. 11301 BM)
- On the fiberscope using the green adaptor included in delivery (Art. no 13 242 XA); On the machine: consult the manufacturer of the washer/disinfector.

If the reprocessing machine cannot be connected to the leakage testing unit for any reason, check for leakage manually as described on page 43 or check the adaptor connection using the 13242 XA.

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For machine reprocessing, attach the green cleaning adapter (Art. no. 11301 CD) for the reprocessing of suction and irrigation channels in cleaning machines to a reprocessing tube and position on the valve connections. Use the cleaning adapter (Art. no. 11301 GG2, see page 8) for reprocessing in ETD2/ETD3 machines. In addition, use the special connecting adapter (Art. no. 11301 KK2, see page 8) for cleaning in ETD3 machines. For fiberscopes with a working channel without a suction valve, connect the tube of the washer and disinfectant up to the working channel ⑧. To avoid pressure losses in the washer and disinfectant, lock unused connections with the LUER-Lock (Art. no. 6011590).

WARNING: Please observe the instruction manual of the washer and disinfectant.

Lay the fiberscope carefully in the tray of the machine. Connect the rinsing tubes to the channels. Proceed with the endoscope cleaning and disinfection program according to the regulations of the machine manufacturer.

CAUTION: After the program has finished, the fiberscopes (including all adapters except 11301 CD) are highly susceptible to corrosion due to the ambient conditions in the reprocessing machine.

Remove the fiberscopes from the machine as soon as possible after the program has finished. Then dry the optical surfaces with a cotton-tipped applicator soaked in 70 % alcohol.

NOTE: If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberscopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.

NOTE: For flexible fiberscopes without channels, a hanging storage position is not essential.

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
Connection for machine reprocessing in ETD3 machines

Dismantling the adaptor 11301 GG2:
Press the button corresponding to each tube and pull the tube off simultaneously.

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
1. Slide the free tube end ① of the adaptor (Art. no. 11301 GG2) onto the olive connector of the green adaptor (Art. no. 11301 CD).

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NOTE: For fiberscopes with a working channel without a suction valve, pull the tube not required, e.g. ① out of the metal block.

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2. For machine reprocessing, attach the green adaptor (Art. no. 11301 CD) to the valve connections of the fiberoscope.

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3. Attach the tube with LUER-lock to the fiberoscope instrument channel ⑧ and secure it.

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
4. Insert the other tube ends in the metal part of the adaptor (Art. no. 11301 GG2).

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5. Connect the adaptor (Art. no. 11301 FF2) to the 'valve for equalizing pressure and connecting leakage tester'.

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
6. Connect the leakage tester to the ETD3 using the adaptor 11301 FF2.

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
7. Connect the cleaning adaptor for cleaning the instrument channel in EDT3 cleaning machines (Art. no. 11301 GG2) to the connector of the supply tube. Then connect up the special connecting adaptor for ETD-3 machines (Art.-Nr. 11301 KK2).

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8. Lay the instrument carefully in the tray of the ETD3 machine.

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9. Connect the rinsing tubes to the ETD3 machine.

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WARNING: To avoid larger pressure losses in the EDT3 machine, connections which are not in use should be locked with LUER locks.

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
10. Screw the adaptor 11301 KK2 onto the designated connections of the ETD3 machine, according to the figure.

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11. Use the cleaning and disinfection program for fiberoptic scopes according to the ETD3 manufacturer's instructions.

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After cleaning has been completed, remove the green adaptor (11301 CD) for cleaning the working channel by pressing together the projections at the side.

NOTE: The green adaptor for cleaning the working channel may be left on the tube for later use.

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Remove the connection from adapter 11301 FF2.

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
CAUTION: The fiberscope could be damaged by liquid which may enter the device via an adaptor for leakage testing (11301 FF2) left on the fiberscope after machine reprocessing. Therefore, you should always dismantle the adaptor (11301 FF2) immediately after cleaning is completed!

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Zerlegen des Adapters 11301 GG2: Press respective button and pull off tube simultaneously.

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
Connection for machine reprocessing in ETD2 machines

Push free end of the adaptor tube (Art. no. 11301 GG2) onto the barbed tube (Art. no. 11301 CD).

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Place the adaptor 11301 CD on the valve nozzle ③ of the fiberscope.

Insert the other end of the tube into the metal part of adaptor 11301 GG2.

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
Place the LUER-Lock on the working channel ⑧ of the fiberscope.

Insert the other end of the tube into the metal part of adaptor 11301 GG2.

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Insert the supply tube with the fitting for ETD2 cleaning machines into the metal part of adaptor 11301 GG2.

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Attach adaptor 11301 FF2 to the valve ⑤ for pressure compensation and the leakage tester connection of the fiberscope.

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Connect leakage tester connection of the ETD2 cleaning machine.

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CAUTION: After the program run, the fiberscopes are greatly at risk of corrosion in the reprocessing machine.

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
The fiberscopes (including all adapters except 11301 CD) must be removed from the machine as soon as possible once the program has finished.

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Then if necessary dry off the optical surfaces with a cotton swab soaked with 70 % alcohol.

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Remove the leakage tester connection from the adaptor 11301 FF2.

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CAUTION: The fiberoscope could be damaged by liquid which may enter the device via an adaptor for leakage testing (11301 FF2) left on the fiberoscope after machine reprocessing. Therefore, you should always dismantle the adaptor (11301 FF2) immediately after cleaning is completed!

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

To disassemble adaptor 11301 GG2: Press respective button and pull off tube simultaneously.

NOTE: If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberoscopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.

NOTE: For flexible fiberoscopes without channels, a hanging storage position is not essential.

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18. Inspection after reprocessing

Before assembly it should be checked that all parts and surfaces are undamaged, free from any organic or chemical residue as well as completely clean and dry.

Replace defective parts or send in for repair. Before being sent in for repair, all parts must be cleaned, wipe-down disinfected and sterilized in the NT procedure if necessary.

If residue is still present after reprocessing, reprocess the parts again, subject them to special treatment, discard or repair. Remove stubborn deposits on glass surfaces with a special cleaning paste (Art. no. 27661).

The part must then be submitted to a complete reprocessing cycle.

Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their new condition, they must be replaced. This applies in particular to seals and rubber caps.

Test the instruments for proper operation: check the telescope for a clear view, the instrument for steerability and all canals for patency.

WARNING: Damage to the tube: The use of ULTRASTOP as an anti-fog solution on the optically distal end faces of the flexible endoscopes can damage the tube connector.

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**19. Gas sterilization/plasma
sterilization**

Sterilization using formaldehyde or ethylene oxide and plasma sterilization are gentle and therefore recommended methods of sterilization.

NOTE: Observe national laws and regulations for the gas sterilization method.

NOTE: Fiberscopes which are to be used for interventions in areas of the body which lack microbial colonies, must be sterilized in sterile packaging and stored in a closed cupboard free from contamination and dust.

CAUTION: If the instrument is to be sterilized, the entire system must be opened by fitting the red cap on the pressure compensation part to ensure pressure compensation.

WARNING: Before commencing gas sterilization, the instrument must be completely dry or else sterile conditions will not be achieved.

- Place the instrument securely on a sterilization sieve.

Ethylene oxide gas sterilization has been validated under the following parameters:

Gas mixture	EO: HCFC*
Temperature (°F, °C)	130 °F ± 5 °F; 54 °C ± 2 °C
Relative air humidity	60 ± 20%
Pressure (p _a , kpa)	58 ... 70 (0.56 ... 0.7 bar/8 ... 10 psi)
Exposure time	120 mins
EO concentration	600 ± 30 mg/L

* Gas mixture consisting of 10 % ethylene oxide and 90 % chlorotetrafluoroethane (HCFC-124; percent by weight)

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For gas sterilization with ethylene oxide, the following airing times must be observed for instruments due to gas absorption of the materials:

	at 20 °C	at 42 °C	at 55 °C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Fiberscopes/ Plastics	5 days	2 days	30 hrs

(minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes.)

Maximum acceptable levels (per 21 CFR 812.100; Federal Register, Vol. 43, No. 122, June 23 1978) of residues following EO sterilization are as follows:

Ethylene oxide:	25 ppm* / 250 ppm**
Ethylene chlorohydrine:	25 ppm* / 250 ppm**
Ethylene glycol:	250 ppm* / 5000 ppm**

** for instruments in contact with skin and mucous membrane,
* for instruments in contact with blood.

The instruments can be used immediately after sterilization with formaldehyde or plasma sterilization. Extra airing time is not necessary.

Remember if you use plasma sterilization that there are limits regarding the size of the lumen and material in regard to instruments which may be sterilized using STERRAD®.

Sterilization using STERRAD NX is also possible.

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Statements as to the efficiency of sterilization for the respective type of instrument (depending on the length and diameter of the lumina) may only be made by the manufacturer.

WARNING: Sterility is only assured if the manufacturer's instructions for the STERRAD® procedure are followed.

WARNING: All deviations from the sterilization parameters recommended for the STERRAD® system must be validated by the user.

CAUTION: Attention must be paid to the information on reprocessing, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

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20. Chemical sterilization

Chemical sterilization (STERIS®) is a special method which is particularly suited for thermolabile goods such as fiberscopes. Except for telescopes with a HAMOU® eyepiece drive mechanism, all KARL STORZ telescopes are suitable for STERIS® sterilization. Color changes to the surface of aluminum alloys could occur during sterilization, which do not however affect the function of the instrument. It must be pointed out in this context that anodized aluminum parts are not always recognizable as such.

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WARNING: Sterility is only guaranteed if the manufacturer's instructions for the STERIS System have been followed. This also applies to use of the sterilization containers.

WARNING: Chemical sterilization is not possible on greased or oiled surfaces.

WARNING: Parts with a lumen must be reprocessed together with the STERIS Quick Connect kit. All parts must be thoroughly cleaned before sterilization using the STERIS PROCESS™, as also applies to the other sterilization methods.

CAUTION: Attention must be paid to the information on cleaning, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

CAUTION: For miniature endoscopes with a removable eyepiece, the eyepiece must either be sealed with a plug or firmly attached to the endoscope before sterilization.

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21. Storage

CAUTION: Do not expose fiberscopes to direct sunlight or excessive heat.

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Sterilized fiberscopes should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage.

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Storage must take place in facilities with normal temperature and air humidity levels.

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You can find special holders for this in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions.

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22. Accessories

Flexible forceps

Carefully wipe the outer surfaces of the flexible forceps with a disposable cloth soaked in disinfectant. Clean the forceps jaws with a brush.

Next, clean in an ultrasound bath. Afterwards, rinse with microbiologically pure/sterile water.

Manual/machine cleaning/disinfection

Both processes can run parallel to decontaminate flexible forceps.

Exception: If the flexible forceps are thermostable, they can be cleaned and disinfected in the 93 °C program.

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For handles which can be dismantled:

Dismantle handle as follows:

- Undo the thumbscrew on the handle and remove the rear part of the handle.
- Undo the thumbscrew on the connecting piece and remove the front part of the handle to the rear.

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Also oil the joints of the jaws by holding the instrument so that the oil can flow into the instrument channel and wet the pulling wire. Then reassemble in reverse order.

The sterilization of flexible forceps is compulsory!

NOTE: Only use special oils whose lubricating effect remains even in autoclaves (art. no. 27656).

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Papillotome/Polypectomy snare

Dismantle the device completely for reprocessing.

- Unscrew thumb handle.
- Undo contact socket.
- Unscrew grip piece.
- Remove inner guide sleeve.
- Pull papillotomy or polypectomy loop out of the plastic probe at the front.

Then reassemble the polypectomy snare in reverse order.

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Both processes can run parallel to decontaminate papillotomy/polypectomy snares.

Exception: If the papillotomy/polypectomy snares are thermoresistant, they can be cleaned and disinfected in the 93 °C program.

Wipe the exterior of the instruments with a cloth soaked in disinfectant. Clean the forceps jaws with a disinfected brush.

Then clean in an ultrasound bath.

Afterwards, rinse with microbiologically pure/sterile water.

Assembly of papillotomes:

- Screw inner guide sleeve into the grip piece
- Fit contact socket.
- Insert loop into plastic probe as follows: insert long end of loop into the rear opening of the probe and push to the rear, keeping the probe stretched.
- Insert shorter end of loop into the front opening of the probe and slide both parts towards the rear keeping them parallel.

In polypectomy snares, the loop end is already shaped. Insert the proximal loop end into the plastic probe from the distal end and push towards rear.

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Flexible injection cannula

For reprocessing, pull the needle out of the flexible casing.

Rinse the cladding, then immerse it in solution and clean in the ultrasound bath.

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Both processes can run parallel to decontaminate papillotome/polypectomy snares.

Exception: If the papillotome/polypectomy snares are thermostable, they can be cleaned and disinfected in the 93 °C program.

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Rinse out the cannula using a syringe with disinfectant solution.

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24. Normatives

24.1 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

24.2 Limitation of liability

As the supplier of this instrument, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions:

assembly, system expansion, readjustment, modification, or repair have been performed by persons authorized by KARL STORZ and that the instrument is used in accordance with the instruction manual.

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The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

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25. Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective instrument. For telescopes, you receive a warranty of 1 year, and for instruments 6 months.

For fiberscopes and equipment, individual repair is necessary. To bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
D-78532 Tuttlingen/Germany

In other countries please contact your local KARL STORZ branch or authorized dealer.

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25.1 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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KARL STORZ Endoscopy China Ltd
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016
People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjil-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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**Important information for users
of KARL STORZ Instruments**

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This manual contains information for personnel authorized to prepare instruments for ORL. The instruments and accessories manufactured by KARL STORZ are, as a rule, sensitive precision instruments which require particular care during preparation. To prevent damage during preparation, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ instruments.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be deduced from the information, figures and descriptions in this manual.

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ABOUT THIS MANUAL.....1

Using these instructions	1
Symbols and hazard information	2
Safety instructions	3
General instructions on cleaning and sterilization	5

OTOLOGY.....12

Suction tubes, aspirator attachment and suction handles	12
Forceps and scissors	12
Hand instruments for the ears	12
Magnetic instrument rack	13
Instrument rack	13
Aural polyp snare	13

INSTRUMENTS FOR RHINOLOGY AND RHINOPLASTY14

Atomizer for the nose	14
Rubber bulb	14
Nasal snare	14
Rasps, saws	16
Flexible tension arm	16

SINOSCOPY, RHINOSCOPY, POSTRHINOSCOPY17

Shaver blade, Paranasal Sinus Shaver Handle, Model 40 711039	17
EC motors 20 711030/32, Drill handles 2514XX	22
STAMMBERGER/SACHSE intranasal drill	23
Micro saws	26
Micro Compass Saw 252100 or Micro Oscillating Saw 252200	27
Punches, rotatable	28
Suction and irrigation sheath	29
Laser rhinoscope	32
LEUNIG/GREVERS laser probe guide instrument	33
Antrum punch, rotatable	34
Biopsy forceps and grasping forceps with suction channel	38
Optical forceps	38

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Adenotomes	39
Mouth gag, Model 744500	39
Mouth gag, Model 744600	40
Salivary duct probe	41
Tonsil snare	41
Laryngeal mirror, standard version	41
BIPOLAR AND MONOPOLAR COAGULATION	42
General information on coagulation instruments	42
Assembly and care	42
Cable connections	43
Forceps for bipolar coagulation	45
Instruments for bipolar coagulation	46
Instruments for monopolar coagulation	48
LARYNGOSCOPY, TRACHEOSCOPY 49	
Distensible oropharyngoscope	49
Distensible laryngoscopes	49
Fiberoptic light carrier, clip-on light deflector	50
Semirigid and rigid, nonseparable forceps, scissors, needle holders	52
High pressure syringe for injecting viscous liquids	53
Tracheoscopes	54
Optical instruments	55
TITAN TELESCOPES AND INSTRUMENTS	56
MR mouth gag	
<i>see Adenotomy, nasopharynx, tonsillectomy, mouth, tracheostomy, larynx section</i>	
Suction tubes	
<i>see Otology section</i>	
Forceps and scissors	
<i>see Otology section</i>	
MAINTENANCE	57
Operating/storage conditions	57
Disposal	57
Directive compliance	57
Limitation of liability	57
Warranty	58
Important information	58
Servicing and repair	58
Repair and replacement program	58

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Using this manual

This manual is not an instruction manual, but instead provides general information regarding cleaning, disinfection, care and sterilization of ORL instrument sets from KARL STORZ.

This manual presumes knowledge of the "Cleaning, sterilization and care of KARL STORZ instruments" manual, in which the various cleaning and sterilization methods are described in detail. If you do not yet have this manual, it can be requested without charge from KARL STORZ.

Instruments whose cleaning, dismantling and reassembly do not require any special knowledge are not listed separately in these instructions. Particularly sensitive product groups, such as telescopes or fiberscopes, have their own manuals. Separate instruction manuals and service manuals are available for units. Other general regulations (Normal use, User qualification, Safety rules when using instruments, etc.) must be heeded.

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Symbols and hazard information

Please read this manual and follow its instructions carefully. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words Warning and Caution stand out more clearly, they are accompanied by a pictogram.

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument.

Note: A Note indicates special information about operating the instrument or clarifies important information.

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Safety information

Warning: The listed cleaning, disinfection and sterilization methods alone do not guarantee disinfection or sterility. This can only be achieved if staff follows recognized and validated preparation methods.

Warning: New instruments are delivered unsterilized and must therefore be thoroughly cleaned, disinfected and/or sterilized prior to initial use and each subsequent reuse.

Warning: For all cleaning of contaminated endoscopes and instruments, the Accident Prevention Regulations (UVV) must be observed (e.g. wearing of protective gloves and safety goggles).

Warning: Remember that sterilization is only then successful if

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Caution: Use of preparation methods other than those listed here may have a negative effect on the function and service life. Any damage resulting from this is not covered by the guarantee.

Caution: The following information on cleaning and sterilization only applies to the thermostable instrument set described in these instructions. For telescopes and fiberscopes, special rules apply (see instructions for telescopes and fiberscopes).

Caution: For the production and application of solutions, precisely observe the information from the manufacturer regarding mixing ratios and immersion times. A maximum immersion time of 60 minutes (including in water) must not be exceeded. Prolonged immersion can lead to damage. After cleaning it is recommended to rinse the instruments with completely demineralized water.

Caution: Constant changing between various preparation methods is especially hard on the materials and must therefore be avoided. A one-off change of method, due to procurement of new cleaning and sterilization equipment for example, is however possible.

Observe the guidelines of the Employer's Liability Insurance Association and equivalent organizations.

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Caution: The instructions for cleaning and sterilization in the manual 'Cleaning, sterilization and care of KARL STORZ instruments' must be observed. The methods for cleaning, disinfection and sterilization are explained there in detail.

Note: Recommended cleaning accessories and sterilization containers can be found in the catalog chapter 'Cleaning, sterilization and care'.

1. Manual preliminary cleaning

The instruments may require manual preliminary cleaning on account of their design (e.g. slits and long, narrow lumens) and depending on the type and degree of soiling.

In addition, a prolonged waiting time before preparation (which allows the contamination to dry on) may make preliminary cleaning necessary while still in the operating room. Organic residues and pharmaceutical products, for example, must be removed mechanically (by brushing, rinsing through with a water gun, etc.), and the instruments must subsequently be rinsed off with completely demineralized or distilled water.

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Tissue remnants sticking or baked to the jaws must be removed with a cleaning brush, and the instrument must be precleaned with a 3% H₂O₂ solution as soon as possible after use. Special care must be taken when cleaning the jaws, as they are precision-engineered parts.

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2. Ultrasound treatment

All non-separable scissors, coagulation instruments and microinstruments must be precleaned with ultrasound (3 - 5 minutes exposure time). The cleaning effect of the ultrasound bath reaches even inaccessible hollow spaces. For this reason, all non-separable scissors and forceps, with or without irrigation channel, rail-sheath instruments, aspirators, circular punches, coagulation instruments and microinstruments must be preliminary cleaned with ultrasound. Ultrasound has no cleaning effect on soft plastic parts because the oscillations are attenuated.

Hollow instruments (aspirators etc.) should be held at an angle in the ultrasound bath so that the hollow spaces can fill. Position large instruments in the bath so that no acoustic shadows are formed. The following items, amongst others, are not suitable for cleaning in an ultrasound bath: all rigid and flexible telescopes, electrical devices, cables, prisms, rubber articles, magnetic discs and fiberoptic light cables.

After cleaning, rinse thoroughly with completely demineralized water and then dry immediately.

Place the instruments in a suitable instrument carrier (e.g. a tray holder).

***Note:** If manual preliminary cleaning is mandatory, the necessary steps are explained in this manual for each instrument*

3. Manual/machine cleaning and disinfection

For preparation, multi-part instruments must be dismantled, and seals, cock plugs and protective caps removed.

Manual cleaning/disinfecting cleaning

For manual cleaning or disinfecting cleaning, use active-cleaning, non-protein-fixing chemicals, where necessary with a verified range of action (DGHM (German Association of Hygiene and Microbiology) listing or comparable).

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Caution: When preparing and using the solutions follow the manufacturer's instructions paying close attention to proper dilution and immersion times.

Note: A list of products which have been tested and shown not to attack the surface and approved by KARL STORZ for manual cleaning and disinfection can be found in the appendix of this manual and on the KARL STORZ home page.

When immersing instruments in cleaning and disinfectant solutions, ensure that air can escape from hollow spaces so that all the internal surfaces are completely wetted. It is best to immerse the instruments at an angle until no more air bubbles can be seen. During cleaning, always ensure that all hollow spaces (sheath lumens, stopcocks, valves etc.) are also cleaned and that there are no obstructions (flow test, for example with a syringe).

Thoroughly clean surfaces, interiors and jaws with a suitable soft brush, sponge or cleaning gun, and rinse. Take care when cleaning small parts.

The cleaning gun with accessories (Art. no. 27660) is ideal for rinsing out hollow spaces with water and drying them with sterile compressed air.

After cleaning, rinse the instruments copiously with completely demineralized water and dry with sterile compressed air.

Thermal cleaning and disinfection in a machine

As far as possible, thermal cleaning/disinfection should be carried out in a machine, followed by steam sterilization at 134°C. The advantage of preparation by machine is that it reduces the risks for personnel and standardizes the cleaning procedure (validatable or validated processes are demanded).

As there is sometimes strong movement of the water, it is necessary to secure all parts in holding devices to avoid damage. To ensure that the instruments are adequately rinsed out, it is essential to connect up the hollow spaces appropriately to the washer/disinfectant.

It is not advisable to use lubricants etc. in the machine because they may cause deposits to form on the instruments and may damage plastics.

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Selective manual care and lubrication of joints etc. is safer and more effective.

It is necessary to perform a final thorough rinse with completely demineralized water in order to remove all chemical residue, as such residue can lead to chemical cross reactions or the formation of deposits during final sterilization.

Note: A list of products tested and approved by KARL STORZ for machine cleaning/disinfection can be found in the appendix of this manual and on the KARL STORZ home page.

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- Success of previous cleaning / signs of residual dirt (as a rule via visual checks of the critical areas)
- Smooth operation of the joints
- General signs of wear
- Exact closing of forceps, clamps and punches etc.
- Secure closing of valves and seals
- Sharpness of cutting instruments
- Condition of telescopes
- Condition of the insulating sheaths of cables, including cold light cables
- Broken fibres in cables

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Prior to sterilization, instruments with a closure mechanism or joint must be treated with lubricants capable of withstanding steam sterilization (e.g. instrument oil Art. no. 27656 B).

Note: Apply a thin coat of special grease to sealing and running surfaces, such as stopcocks and plungers (Art. no. 27657). This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the stopcocks leakproof. Oil is not suitable as a sealant.

Warning: Only use special white oil-based lubricants as the sterilant is not able to penetrate all types of oil in order to sterilize the underlying surfaces.

Note: Leave the stopcocks dismantled for sterilization and only insert them into the housing once sterilization is complete.

5. Sterilization

For almost all instruments, steam sterilization by the fractionated vacuum method at 134°C (3.04 bar.), gas sterilization (formaldehyde, ethylene oxide), plasma sterilization and chemical sterilization (STERIS) are possible. Steam sterilization at 134°C is usually the gentlest and most economical method of sterilization. This also applies to rigid telescopes because the longer exposure time for steam sterilization at 121°C can have a negative effect.

Sterilization is only ever possible on surfaces that have been cleaned.

Almost all instruments can be sterilized without being dismantled. Only tighten screws and closures gently, so that they are easy to undo later.

Note: If the design of the instrument necessitates sterilization in the dismantled state to ensure effectiveness, this is expressly stated for the article in question in this manual.

For steam sterilization (with saturated steam at 134°C + 3°C), the process defined in EN 285 or DIN 58946 Part 1, Number 3.25.2 should be used.

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Use prevacuum methods for instruments with a lumen. The steam must fulfill the requirements of EN 285, Appendix B, in order to avoid discoloration and stains.

These procedures are also permissible for telescopes labeled 'AUTOCLAVE'.

Caution: When sterilizing telescopes, please follow the instructions in the Telescopes manual.

Warning: Long lumens, slits and channels, and assembled instruments can only be sterilized without restrictions with the fractionated prevacuum method.

Note: If it is necessary to deactivate the pathogens responsible for Creutzfeldt-Jakob disease, the steam sterilization method recommended by the WHO can be used (18 min at 134°C), the sterilization time may be extended to up to an hour at 134°C in certain cases.

As there are various types of cleaning, disinfection and sterilization, it is essential to ensure that recognized and validated (tested) preparation methods are used.

Refer to the manual 'Cleaning, sterilization and care of KARL STORZ instruments' for more detailed information.

6. Storage/Transportation

To avoid damage during transportation and storage, place the instruments carefully in suitable holding devices.

Suitable accessories, storage systems, and sterilization trays for special instrument sets can be found in the relevant chapter of every KARL STORZ catalog.

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Suction tubes, aspirator attachment and suction handles


Cleaning

Clean the exterior with a sponge or brush. Thoroughly clean the lumen with a brush or wire stylet, then rinse thoroughly with water and dry with sterile compressed air.

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Forceps and scissors


Cleaning

Carefully clean the jaws with a soft brush. With rail-sheath instruments, open the flap over the sheath and release the slide. Then clean with a soft brush whilst open. Rinse with water and dry with sterile compressed air. Clean in an ultrasound bath.

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
Oil the joints (Instrument oil 27656 B). Whilst doing this the instrument is to be held in such a way that the oil is able to flow into the sheath and also cover the wire pull.

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Hand instruments for the ears

The hand instruments for the ears need to be stored in a particular way to prevent the instruments tips from being damaged (e.g. instrument rack, magnetic instrument rack, knobbled silicone mat).


Caution: Magnetic instrument racks must not be cleaned in ultrasound baths.

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Magnetic instrument rack


! Caution: Magnetic instrument racks must not be cleaned in ultrasound baths.

Both machine cleaning and steam sterilization are possible. Ensure that the instrument is not subject to vibrations when stored in order to prevent abrasion damage to the magnets.

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Instrument rack

Instrument racks may be cleaned in ultrasound baths. Following machine cleaning and steam sterilization, it must be ensured that the instrument holders in the instrument rack are dry, and that the joints of the instrument rack are treated with steam sterilizable silicone oil prior to steam sterilization. Any liquid residue is to be dried off with sterile compressed air.

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
Instruments for rhinology and rhinoplasty

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Atomizer for the nose

Remove the cap with atomizer tube from the glass bulb. Then rinse the glass bulb and channels with completely demineralized water and dry with compressed air. Rinse the disinfected atomizer with sterile water thoroughly and dry with compressed air. Then steam sterilize the dismantled atomizer at 134 °C.

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
Rubber bulb

Only the exterior of the rubber bulb may be cleaned; no water must penetrate into the bulb during cleaning. Steam sterilization at 121°C is gentler than sterilization at 134°C, it nevertheless results in a significant amount of wear.

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Nasal snare


Disassembly

Cut through the snare wire using suitable scissors.

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Unscrew screws ① and ② (do not remove completely), and remove the guide tube from the finger grip. Pull out the loop guide.

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Cleaning

Thoroughly brush and rinse the guide tube and loop guide. Dry with sterile compressed air.

Disinfection


Immerse the dismantled and cleaned nasal snare in disinfectant solution. Then rinse thoroughly with completely demineralized water and dry with sterile compressed air.

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Assembly


Apply a few drop of instrument oil (Art No. 27656 B) to the loop guide and insert into the guide tube from proximal. Place the guide tube on the finger grip. The projection on the guide tube must fit in the recess of the finger grip.

Pull the loop guide out so that it reaches into the opening below screw ②.

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Tighten screw ①, ensure that the loop guide is in a horizontal position when doing this.

Tighten screw ② on the finger grip. Bring the finger grip into the front position.

Ensure that the distal openings of the guide tube are congruent with the distal openings on the loop guide.

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
Insert the snare wires into the two front openings of the loop guide from right and left. The loop guide is inserted into the guide tube by pulling back the finger grip.

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Rasps, saws


Caution: Manual cleaning may only be performed with plastic brushes. The cutting edges may be damaged if metal brushes are used.

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Flexible tension arm

Remove the rotating clamp ① and adaptor ② if necessary, and relieve the load on the arm.


Cleaning is performed in an ultrasound bath. Then rinse with completely demineralized water, disinfect, and sterilize in an autoclave.

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Irrigation

If the shaver is given to nursing staff or laid down, the attachment should be briefly rinsed out with saline solution with the cutting opening open, in order to prevent the attachment from becoming blocked.

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Disassembly

Pull the shaver attachment out of the instrument holder of the shaver handle and remove the suction tube from the shaver handle.

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
Model 407110 39:

Disconnect the shaver handle housing (1) from the shaver handle motor (2) by pressing the two locking buttons (3).

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Dismantling the shaver attachment

The inner part of the shaver attachment can be pulled out toward distal.
In the case of curved attachments, the inner part is flexible.

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Caution: Do not mix up the individual inner and outer shaver attachments. If incorrectly assembled, friction at the tip of the shaver attachments can lead to an intense build-up of heat.

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
Cleaning/sterilization

Firstly, clean the shaver handle and the used attachments mechanically with water or a cleaning and disinfectant solution, and a soft brush. Do not use hard or sharp cleaning materials! After brushing, rinse out the handpiece suction channel and the attachments with a strong jet of water until no more residue can be seen.

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
In order to avoid material damage, carry out the final rinse using completely demineralized water. Then dry using sterile compressed air. The handle and blades can now be gas sterilized or steam-sterilized at 134°C (+3°C).

Note: After manual preliminary cleaning, if necessary, we also recommend cleaning and disinfecting the shaver handle and blades by machine. In the case of the attachments, however, you must ensure to rinse the sheaths (by connecting them specially in washers and disinfectors). The shaver handle housing can be connected via a LUER lock connector d to washers and disinfectors.

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
* **Before** each sterilization the shaver handle housing must be sprayed with universal spray 280052 B at the three points marked on the handle.

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Assembly

Connect the shaver handle housing (1) to the shaver handle motor (2) via the socket pins. Connect the suction tube to the shaver handle.

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
Insert the inner sheath of the shaver attachment into the outer sheath.

Note: To avoid damage to curved shaver attachments, the cutting opening of the inner sheath must be level with the orientation pin of the set wheel.

Note: Shaver attachments must be lubricated prior to sterilization (oil 27656 B).

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Press the shaver attachment into the instrument holder of the shaver handpiece.

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Instructions for use

With model 407110 40, the set wheel has to be fixed manually. Always leave the tip of your index finger in one of the notches on the set wheel. This allows the cutting direction to be set and fixed at any time.

Warning: If the set wheel is not fixed, the patient may be injured by the uncontrolled rotation of the shaver attachment.

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Blockage of the shaver attachment

In the event of a possible blockage in the attachment, remove this and clean/rinse with a syringe.

Note: In order not to impair the performance of the shaver system, ensure that there is constant suction, preferably via a standard surgical aspirator.

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EC motors 20 7110 32 Drill handles 2524XX

Dismantle the EC motor from handpiece and burrs for cleaning.

The EC motor can be cleaned in the washing machine and steam sterilized at 134°C or sterilized with gas.

Note: Rinse out the cooling/irrigation tube with distilled water after every use.

For cleaning, the chuck must be open. If necessary, turn the clamping collar to the 'Open' position. The handpiece and burr are cleaned mechanically with warm water or a cleaning and disinfectant solution (see Appendix) and a soft brush. Do not use hard or sharp cleaning materials. In order to avoid material damage, carry out the final rinsing with completely demineralized water. Finally dry with sterile compressed air.

The individual parts can be steam sterilized at 134°C, sterilized with gas or prepared in a washer/disinfector.

Before every sterilization the handpiece and the EC motor must be treated with universal spray 280052 B, or at least twice a day (if it is used frequently without sterilizing in between).

Note: The clamping collar (chuck) must be open.

To do this, place the handpiece/EC motor on the spray nozzle and lightly press it onto the spray nozzle. Cover the handpiece/EC motor with a soft cloth or cotton wadding and actuate the spray for a period of approximately one second.

Caution: Do not kink the motor cable.

Warning: Let the EC motor cool down after autoclaving.

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Dismantling

Unlock and remove the outer sheath.
Unlock and remove the guide sheath.

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

Press the release button (5) on the sliding handle downward and pull the grinder off toward distal.

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Note: There are two release buttons on the sliding handle (one on each side). Since they are connected (rocker operation), it doesn't matter which one is actuated.

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Dismantle the valve on guide sheath 660000 T. To do this, release the thumbscrew and remove the valve.

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
Cleaning//Disinfection

Clean the exterior of the instrument with a soft sponge or soft brush. Rinse with completely demineralized water and dry with compressed air.

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Carefully clean the valves and connectors with a brush. The valve connectors and sheath lumen must then be thoroughly rinsed with completely demineralized water and subsequently dried with compressed air. Clean/disinfect the cleaned guide sheath as well as the disassembled valve parts in an ultrasound bath. Finally, rinse the cleaned and disinfected parts with sterile completely demineralized water and dry with sterile compressed air.

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Assembly

Lightly grease the running surfaces of the valve (special grease for valves: Art. no. 27657).

Note: This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the valves leak-proof. Oil is not a suitable sealant.


Insert the valves and fix with the spring cap. Dismantle again for sterilization.

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Insert the grinder. Ensure that the wide-surfaced part of the grinder is facing downward.

Push the grinder backward until it locks into the opening on the handle part next to the release button ⑤.

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
Attach the guide sheath and lock into place. (observe the '0' marking).

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Micro saws**Cleaning**

The micro saw is to be separated from the EC motor for cleaning.

The EC motor may be steam sterilized at 134°C or gas sterilized.

Remove the union nut and saw blade.

Mechanically clean the micro saws and saw blade with warm water or a cleaning and disinfectant solution (see appendix) and a soft brush. Do not use hard or sharp cleaning materials! In order to avoid damaging the material, use completely demineralized water for the final rinse. Dry with sterile compressed air.

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**Micro Compass Saw 252100 or
Micro Oscillating Saw 252200****Cleaning**

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The EC motor can be steam sterilized at 134°C or gas sterilized.

Remove the union nut and saw blade.

Mechanically clean the Micro Compass Saw 252100 or Micro Oscillating Saw 252200 with warm water or a cleaning and disinfectant solution (see appendix page 59) and a soft brush. Do not use hard or sharp cleaning materials! In order to avoid damaging the material, use completely demineralized water for the final rinse. Dry with sterile compressed air.

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Sterilization/Disinfection

The individual parts can be steam sterilized at 134°C, gas sterilized or immersed in disinfectant solution or prepared by machine.

Before every sterilization, the micro saw and EC motor must be treated with Universal Spray (see accessories).

To do this, place the micro saw/EC motor on the spray nozzle and lightly press it onto the spray nozzle. Cover the micro saw/EC motor with a cloth or cotton wadding and actuate the spray for approx. 1 second (max. 3 seconds).

Following steam sterilization, the vacuum-sealed micro saw must be dried in a drier or allowed to dry for 24 hours at room temperature.

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
Punches, rotatable

Disassembly

Release the thumbscrew on the handle part by turning clockwise, and unscrew the working insert.

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
Cleaning

Clean the disassembled parts with a brush in an appropriate cleaning solution and then rinse well with completely demineralized water. Dry with sterile compressed air.

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Disinfection

Immerse the dismantled and cleaned punch in disinfectant solution. Then rinse well with sterile completely demineralized water and dry with sterile compressed air.

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Assembly

Insert the slide in the seat of the working insert and push forward (distal).

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
Push the handle together, ensuring, however, that a small gap remains (see Fig. 1, 'S' marking). Screw the working insert onto the handle in a counterclockwise direction. The working insert can be fixed into place by twisting the thumbscrew in a clockwise direction.

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- ① Suction/irrigation sheath, can be rotated
- ② Telescope channel cleaning adaptor (Art. No. 723540 L)
- ③ Press Button Valve, bistable: (Art. No. 723620 DV) Actuation switches irrigation on/off completely.
- ④ Push Button Valve, monostable: (Art. No. 723621 SV) Irrigation only as long as the valve is actuated.
- ⑤ Handle
- ⑥ Suction port
- ⑦ Irrigation port
- ⑧ Irrigation channel cleaning adaptor (Art. No. 723540 K)
- ⑨ Sealing rings (replacement seals for valves, set of 2x5: Art. No. 7249291)

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Disassembly


Screw out the valve from the handle ⑤.

Check all the sealing rings ⑨.
Replace any damaged sealing rings, observing the different sizes of sealing ring.
Replacement seals for valves:
Set of 2x5; Art. No. 7249291

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Unlock and remove the telescope.
Pull the suction/irrigation sheath out of the handle ⑤.

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
The cleaning adaptors are suitable for use in the cleaning machine or for rinsing out manually with a syringe.
Preliminary cleaning to be carried out immediately after use: With active irrigation rinse for at least 10 seconds. This prevents impurities from drying out.
Clean the suction and irrigation sheath ① and handle ⑤, if necessary, in the ultrasound bath.

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Thoroughly clean the suction and irrigation sheath lumen ① with a brush. Carefully clean the surfaces with a sponge or soft brush.
Put and fix the cleaning adaptor ② in the telescope channel. Push the distal side of the cleaning adaptor ⑧ onto the suction and irrigation sheath ① as far as it will go. Connect to the cleaning machine or rinse out manually with a syringe.

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
Press button valve: clean with button pressed (in the irrigation position).
Push button valve: keep button pressed during irrigation.

Then rinse all the parts thoroughly with completely demineralized water and dry with sterile compressed air.

Press button valve: apply a drop of instrument oil (Art. No. 27656 B) to the opening.

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Assembly


Ensure that all the sealing rings ⑨ are fitted to the valve.

Insert the valve in the handle ⑤ and tighten up lightly.

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
Insert the suction and irrigation sheath ① into the handle ⑤. The suction and irrigation sheath can be rotated in the handle.

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Insert and lock the suitable telescope.
Connect up the suction and irrigation tubes.

Warning: Be sure never to confuse the suction and irrigation ports!

Warning: Observe the instruction manuals for the irrigation and suction unit!
Check the instrument to make sure it is in good working order.

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
Laser rhinoscope


Remove the telescope. Release the sheath locking mechanism and remove the laser application insert from the sheath.

Remove the rubber caps from the working channel. Release the thumbscrew.

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Cleaning

Clean the sheath with a sponge or brush. Thoroughly brush the telescope channel and working channel.

Disinfection


Thoroughly rinse the sheath and the telescope and working channel with completely demineralized water. Next clean/disinfect the instrument in an ultrasound bath. Rinse well with sterile completely demineralized water and dry with sterile compressed air.

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Release the thumbscrew and remove the valve on the rhinoscope sheath.

Carefully clean the valves and connectors with a brush. The inside of the connectors and sheath lumen must be rinsed thoroughly each time they are cleaned, and then dried with sterile compressed air.

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
Reassemble the valves of the laser rhinoscope sheath. Insert the control instrument into the sheath and lock into place.

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**LEUNIG/GREVERS laser probe guide
instrument**

Remove, and if necessary replace, the rubber cap from the probe channel.

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
Cleaning

Thoroughly brush and rinse the suction channel, dry with sterile compressed air.

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
Apply a few drops of instrument oil (Art. no. 27656 B) to the distal end of the instrument.

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
Disassembly

Rotate the knurled nut on the forceps tip in the direction of the arrow and hold firmly in the end position.

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

Remove the punch attachment by holding at a slight angle and moving in the direction of the arrow (do not lever out).

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Assembly of the cleaning adaptor on the punch attachment

Push the cleaning adaptor 459015 LL onto the punch attachment; the attachment clicks automatically into place. Ensure that the punch attachment projection is lined up with the cleaning adaptor slit.

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
Cleaning

Clean the exterior surfaces of the handle and punch attachment with a sponge or brush.

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
Carefully clean the jaws of the punch attachment with a brush.

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
Attach the cleaning adaptor to the punch attachment.

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
Rinse out the punch attachment with a cleaning gun.
Dry with sterile compressed air.
Remove the cleaning adaptor again.

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
Oil the joint of the punch jaws. Whilst doing this the instrument is to be held in such a way that the oil is able to flow into the sheath and also cover the wire pull.

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

Assembly

Open the handle as far as it will go.
The slit for mounting the punch attachment must
point upward toward 12 o' clock.

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
Attach the handle to the punch attachment at
6 o'clock, as shown in the picture. When doing
this, it is essential to ensure that the punch attach-
ment projection points upward toward
12 o'clock.

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
The handle and punch attachment are to be
aligned; the punch attachment is automatically
fixed into place when the handle is closed.
Finally fix the thumbscrew.

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
**Biopsy forceps and grasping forceps
with suction channel**

Clean the exterior with a sponge or brush. Carefully
free the jaws from any tissue residue. Thoroughly
brush the suction channel, rinse with water and
blow dry using sterile compressed air. Oil the mov-
able parts.

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
Optical forceps

Clean the exterior with a sponge or brush.
Carefully free the jaws from any tissue residue.
Thoroughly brush the telescope channel, rinse with
water and blow dry using sterile compressed air.

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
Oil the jointed parts of the forceps jaws. Whilst
doing this the instrument is to be held in such a
way that the oil is able to flow into the joint and
also cover the wire pull.

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The joints on the handle part must also be oiled.

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Adenoid curette

Clean the instrument with a plastic brush in an appropriate cleaning solution, and then rinse well with completely demineralized water. Dry with compressed air.

Inspect the cutting edge for damage, and, if necessary, repair.

Note: Repairs may only be carried out by KARL STORZ, as the edge of the product requires special handling.

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Adenotomes

Remove the lower part of the adenotome from the main section.

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
Release the blade by lifting it out of the anchor mechanism and carefully pull it out of the front part of the adenotome. Inspect the blade for damage and, if necessary, replace.

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Mouth gag, Model 744500

Remove the silicone tube from the tooth hook and clean with a brush. Inspect for damage, and, if necessary, replace.

Clean the instrument with a plastic brush in an appropriate cleaning solution, and then rinse well with completely demineralized water. Dry with compressed air.

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
Apply a drop of instrument oil to the catch spring.

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Mouth gag, Model 744600


Release the thumbscrew and pull out the wire handle. Remove the silicone rubber from the wire handle and clean with a brush. Inspect for damage, and, if necessary, replace.

Clean the instrument with a plastic brush in an appropriate cleaning solution, and then rinse well with completely demineralized water. Dry with compressed air.

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Applying some silicone oil to the wire handle makes it easier to pull on the silicone tube.

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
Salivary duct probe


Note: The probe is made of silver, it is therefore possible that the surface may become discolored due to oxidation.

Discoloration can be removed using DURAGLIT polishing wool (Art. no. 27658).

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
Tonsil snare

Release the catch on the ring handle by pressing together. Keep the catch pressed and push the catch and ring handle all the way forward. Remove the tonsil snare.

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
Release the knurled screw and remove the snare attachment. Clean with a brush and rinse well with water. Dry with sterile compressed air.

To assemble, proceed in the reverse order.

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Laryngeal mirror, standard version

Caution: The mirror mounts are not leakproof. Only clean with a damp cloth or sponge. Do not immerse in cleaning, disinfectant or irrigation solution. Only steam sterilize in autoclaves at 134°C.

Note: The more the mirror is used, the quicker it will become matt (blind), and must then be replaced (This does not apply to the rhodium mirror).

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Bipolar and monopolar coagulation

General information on coagulation instruments

Cleaning of coagulation instruments during surgery

The coagulation surfaces on coagulation instruments should be wiped off with a moist compress between uses during surgery. Incrustations may be loosened by briefly immersing in a 3% H_2O_2 solution. Before reusing, rinse with sterile water and wipe off again.

Manual cleaning

Remove sticky tissue remnants on the tip carefully with a cleaning brush.

Sheath tubes with lumen must be rinsed out well on the inside.

Dry with sterile compressed air, particularly in the area of the plug.

Note: The HF cable and electrode connectors must be clean and dry. Moisture and dirt or residue from disinfectants may lead to malfunctions during operation.

Assembly and care

Incidents during high-frequency surgery which have resulted in damage to the instrument set or even injuries to the user or patient, may be traced back, under closer examination, to incorrect use or defective equipment. A frequent source of errors here proved to be damaged cable and/or electrodes. A further cause may lie in the use of unsuitable accessories from other manufacturers which are not compatible with KARL STORZ units. Safe and proper use of KARL STORZ instruments is only possible using KARL STORZ high-frequency cables and electrodes.

When testing for proper operation, check the insulation material especially for damages such as cracks and chipped material.

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Cable connections

Warning: The use of damaged HF cables may lead to very serious dangers. Check the cable before every use for visible damage. Damaged cable must never be reused.

HF cables are subject to a certain degree of wear, particularly with improper handling and maintenance. For example, damage to the cable terminal or cable insulation may lead to HF voltage flashovers on the working element and therefore to the metallically conductive parts of the instrument. The parts of the instrument at HF potential then become dangerous during patient application to the patient and surgeon.

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- burn marks on the cable terminal
- defects on the connectors
- breaks in the cable or cable insulation
- missing pin for closing the second hole on the cable terminal.

If damage is found or suspected on HF cable, it must be discarded.

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Cleaning high frequency cables Preparation by machine

HF cables can be prepared by machine both chemically and thermally using special programs for sets of endoscopic instruments.

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Note: We recommend using completely demineralized water at all stages of the program. The selection of chemicals and their ratios must also be observed precisely. The pH values should be in the neutral range. Protective agents should not be used in the machine because coatings may form on the metal contacts and they may damage plastics.

Caution: After immersing the cable in solution carefully dry it out in order to avoid any malfunctions due to residual moisture at the mating contacts. After cleaning, steam sterilization should preferably be performed at 134 °C by the fractionated prevacuum method.

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Sterilization

The HF electrodes can be sterilized using steam at up to 134°C ± 3° and, if necessary, for up to 18 minutes. They may also be sterilized with gas or plasma. Follow the sterilizer manufacturer's instructions exactly.

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Disinfection

Clean/disinfect in an ultrasound bath. Then rinse with completely demineralized water and dry with sterile compressed air.

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Caution: To prevent damage to the insulation layer, the instrument must not be placed loose with the other instruments in the cleaning tray, but fixed into place.

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Always use KARL STORZ HF connecting cables which are specially designed for use with forceps for bipolar coagulation:

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
Instruments for bipolar coagulation

Disassembly

Twist the electrode connection socket ① through a ¼ turn in a counterclockwise direction, and then carefully remove the electrode.

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Manual cleaning

Remove any tissue residue which is attached to the electrode with a soft cleaning brush. Clean the interior of the sheath tubes with lumen with a brush or stylet and rinse thoroughly. Then, rinse thoroughly with completely demineralized water and dry with sterile compressed air, paying particular attention to the contact connections.

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
Disinfection

Clean/disinfect the suction channel in an ultrasound bath. Do not immerse the electrode in the ultrasound bath, rather clean/disinfect with a soft cloth. Then rinse all parts with completely demineralized water and dry with sterile compressed air.

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Assembly

Carefully insert the electrode in the suction tube and fix electrode connection socket ① into place by twisting through a ¼ turn in a clockwise direction.

Caution: The instruments for bipolar coagulation must not be sterilized with hot air.

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Always use KARL STORZ HF connecting cables which are specially designed for use with instruments for bipolar coagulation:

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26002 M, 26004 M, 26005 M, 26006 M,

Instruments for monopolar coagulation

Cleaning

During surgery, wipe down the coagulation surfaces of the instrument between uses, using a moist surgical pad. Remove any incrustations by briefly immersing in a 3% solution of H₂O₂. Then rinse with sterile water and wipe down again.

Manual cleaning

Remove any residues of tissue clinging to the electrode with a soft cleaning brush. Sheath tubes with a lumen must be rinsed out thoroughly. Finally, dry with compressed air, particularly the contacts.

Disinfection

Clean/disinfect in an ultrasound bath. Then rinse thoroughly with completely demineralized water and dry with sterile compressed air.

Caution: The forceps for bipolar coagulation must not be sterilized with hot air.

Note: The connectors on the HF cable and electrodes must be clean and dry. Any moisture, dirt or residues of disinfectant may cause malfunctioning during operation.

Always use KARL STORZ HF connecting cables which are specially designed for use with instruments for monopolar coagulation:

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**Distensible surgical laryngoscope
and oropharyngoscope**

Carefully clean the exterior surfaces with a sponge or brush.

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Clean the lumen with a flexible brush. Rinse with
completely demineralized water, dry with sterile
compressed air.

When machine cleaning, it must be ensured in
particular that the inner areas are rinsed properly
(e.g. by attaching onto the nozzle in the MIS trol-
ley).

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
The oropharyngoscope can be cleaned/disinfected
in an ultrasound bath. Then rinse with completely
demineralized water, dry with sterile compressed
air.

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
Distensible laryngoscopes

Disassembly

Undo knurled screw ① and remove light carrier.

On instruments with telescope holder: undo knurled nut ② and remove telescope holder.

Caution: Hold with the handle pointing upward because the knurled nut contains individual parts.

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Cleaning


Clean the external surface with a sponge or brush. Rinse with water. Dry with sterile compressed air.

Caution: The light carrier must not be cleaned in an ultrasound bath.

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Assembly

Attach light carrier. Oil moving parts. If telescope holder is fitted: Position ball of telescope holder ③.

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
The individual parts of the telescope holder are assembled in the following order (from bottom up):

- Ball ③
- Brake cylinder, comprising
 - sleeve ④
 - rubber bolt ⑤
- Sleeve ⑥
- Screw cap ⑦

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**Fiberoptic light carrier, clip-on light
deflector**

Cleaning

Remove the light carrier connection adaptor ①. Completely immerse in cleaning/disinfectant solution, and allow to take effect. Clean with a sponge. Rinse with completely demineralized water and dry with sterile compressed air. Clean the light entrance and exit faces with a soft cloth, sponge or a cotton tipped applicator soaked in 70% alcohol.

Mount the light carrier connection adaptor ①.

Note: The fiber areas in the light inlet piece in particular must be carefully dried with alcohol. Any disinfectant and cleaning agent residue in the light inlet piece can become burnt on when the light carrier is closed, this can seriously impair light transmission.

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Machine cleaning and disinfection

The choice of rinsing program must be agreed with the machine manufacturer. In most cases special programs for endoscopic instrument sets (sometimes also for telescopes) are available which are more suited for cleaning and which are gentler on materials than general programs.

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
Sterilization

Sterilize the cleaned and disinfected instrument in an autoclave. Use suitable sterilization containers for optical instruments. On completion of the sterilization cycle, the fiberoptic light carrier must be allowed to cool down to room temperature in the sterilization container before the container is opened.

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Semirigid and rigid, nonseparable forceps, scissors, needle holders

Cleaning

Clean the sheath of semirigid and rigid instruments with a sponge.


In the case of instruments with a suction channel, the channel must be carefully rinsed out.

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Carefully clean the jaws with a brush. Next, clean/disinfect in an ultrasound bath. Then rinse off with completely demineralized water and blow dry with sterile compressed air.

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
Finally oil the joint of the forceps jaws. Whilst doing this the instrument is to be held in such a way that the oil is able to flow into the cladding and also cover the wire pull.

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Also oil the movable parts of the handle.

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High pressure syringe for injecting viscous liquids

Disassembly

Remove the injection cannula 27200 S/T from the plastic syringe. Remove the plastic syringe from the retaining sleeve and dispose (disposables)

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Cleaning

Clean all parts (with the exception of the plastic syringe) with a soft sponge or a soft brush in appropriate cleaning solution. Clean the injection cannula, if necessary use a suitable stylet for this purpose.

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Then rinse out the parts thoroughly with completely demineralized water and dry with sterile compressed air.

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
Apply a drop of instrument oil to the movable parts.

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
Cleaning

Clean the accessories with a soft brush. Remove any sealing caps ①, and replace if necessary. Thoroughly brush the channels, rinse with water and dry with sterile compressed air.

Clean the exterior with a sponge. Brush the lumen thoroughly. Rinse out with completely demineralized water and dry with sterile compressed air.

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
Clean/disinfect all parts, with the exception of the sealing caps, in an ultrasound bath. Then rinse with sterile completely demineralized water and dry with sterile compressed air.

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Optical instruments

Cleaning/Disinfection

Clean the exterior of the forceps with a sponge, and the jaws with a brush.

Thoroughly brush the telescope channel. Then clean/disinfect in an ultrasound bath.

Rinse with water and dry with sterile compressed air.

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
Oil the jointed parts of the forceps jaws. Whilst doing this the instrument is to be held in such a way that the oil is able to flow from the distal end into the joint and cover the wire pull.

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Operating/storage conditions

	Temperature	Rel. humidity
Storage	-10°C ... +60°C	10% ... 90%
Operation	+10°C ... +40°C	30% ... 70%

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Disposal

Standard batteries and rechargeable battery inserts may contain toxic substances so dispose of in accordance with the environmental regulations.

To dispose of the instrument set, no special measures are necessary.
National laws and regulations must be observed.

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Directive compliance

The medical products described in these instructions bear the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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Limitation of liability

As the supplier of these medical products we only regard ourselves responsible for the safety, reliability and performance of the products if assembly, system expansions, readjustments, modifications, or repairs have been performed by KARL STORZ, and the product has been used in accordance with its operating instructions at all times.

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Warranty

For **two** years as of delivery to the end customer, we agree to replace the goods free of charge if proof can be provided of faulty materials or faulty workmanship. In doing so we cannot accept the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
D-78503 Tuttlingen/Germany

Any repairs or modifications to medical products by unauthorized persons shall release us from any liability for their safety. Any such repair or modification performed during the warranty period shall void all warranty.

KARL STORZ gives no warranty on the correct functioning of units or instruments which have been repaired by unauthorized third parties.

Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and units must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/units to the sender.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine KARL STORZ parts.

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For telescopes, you receive a warranty of 1 year, and for instruments 6 months.

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For fiberscopes and equipment, individual repair is necessary.

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Typ

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths

Note: You will also find the current list of approved products on our Web site (www.karlstorz.de or www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Chemicals for the preparation of telescopes and endoscopic instrument sets
Manual cleaning and disinfection/ultrasound

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU telescopes with ocular focusing, magnifiers with focusing, motors and their handles, IMPERATOR drill handles, prismatic light deflectors. The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to proper dilution and recommended immersion times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Endo GmbH & Co. OHG	EndoStar	1,2,3
	EndoStar forte	1,2,3
Ester Pharma GmbH	EndoStar	1,2,3
	EndoStar forte	1,2,3
Fresenius AG	MediStar	1,2
	MediStar forte	1,2
Hofla Pharma AG, 20.0	Palbert Hebla	1,2
Indaba Industrie E Commerce LTDA	Lebalesto	2
Jose Collado S.A.	Dardor 4000 Liquido	1
	Dardor 9000	2
	Dardor Sinalehd 2000	2
	Aluox	1,2
Lysolform Dr. Rosemann GmbH	Desclon	1,2,3
	Lysolform 3000	1,2
MediStar International	MediStar	2
	MediStar forte	1
	Edisonite Super	1,3
Merz Hygiene GmbH	Macadent-4S	2,3
	Macadent-Zumskov	1,3
	Macocot-T neu	1,2,3
METASYS Medizintechnik GmbH (formerly LONZA)	ID 50	2
Orochemie	A 20	1,3 A
Promocent AB	Wavacite	2
	Desclon Extra	1,2
Schumacher, Dr. GmbH	Endo	2
	Endo forte	2,3
	Endo N	1,2
	Endo TB	1,2,3 A
	Groenep AF forte	2,3
Schülke & Mayr GmbH	Groenep Instru AF	1,2,3 A
	Groenep Med forte	1,2
	Groenep	1,3
	Lytefol V	2
Starm	EndoCare 2	1,3
	neodisher LIM 2	1,3
	neodisher mediclean	1
	neodisher mediclean forte	1,3
Wegert, Dr. GmbH & Co.	neodisher medizim	1,3
	neodisher SeptoClean	1 A
	neodisher SeptoPreClean	1,2
Whitney Industries PTY Ltd	Aldal Plus	2
	Matra	1

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set.

For **flexible endoscopes**, only immerse the examination sheath, **not** the housing/handle part.

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The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU telescopes with ocular focusing, magnifiers with focusing, motors and their handles, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
Steris	Thermosept RKN-zym	1
	Instru-Klenz	1, 3
	neodisher DuoClean	1, A
	neodisher FA	1
Weigert, Dr. GmbH & Co.	neodisher FA forte	1, 3
	neodisher Oxivarlo	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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To be filled out by instrument owner:

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KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSKOPIY Kasachstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@it.te.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 – 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

Karl Storz Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
#06-07 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

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KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay, Hong Kong,
People's Republic of China
Phone: +8 52 28 65 2411
Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725
Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188
Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175
Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281
Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6F, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966
Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

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**Important information
for users of KARL STORZ
instruments**

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These instructions contain information for personnel authorized to prepare instruments for gynecology. The instruments and accessories manufactured by KARL STORZ are, as a rule, sensitive precision instruments which require particular care during preparation. To prevent damage during preparation, support is provided here for disassembly, cleaning, disinfection, assembly and sterilization of the KARL STORZ endoscopic instruments.

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KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations and descriptions in these instructions.

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**Important information for users of
KARL STORZ instrumentsU1**

About this manual	4
Using this manual	4
Symbols and hazard information	5
Safety instructions	6
General instructions on preparation of medical devices.....	10

**Hysteroscopes for examination and
surgery..... 17**

Sheaths with LUER lock stopcock.....	17
Semirigid instruments	19
Flexible instruments	20
BETTOCCHI polypectomy snare	22
Compact miniature hysteroscopes	24
BETTOCCHI B.I.O.H Compact Hysteroscope	26
CAMPO compact instrument set.....	30
Overview of preparation.....	33
HF loops and electrodes, see chapter on <i>'Resec- toscopes for intrauterine HF electrosurgery'</i> Flexible hysteroscopes, see separate <i>'Fiberscopes' manual</i>	

**Resectoscopes for intrauterine HF
electrosurgery.....34**

Safety instructions for the use of resectoscopes.....	34
Resectoscopes model 26040xx / 26050xx / 26055 xx	35
Resectoscopes model 26040 xx / 26040 xx / 26055 xx.....	37
HF electrodes: Safety instructions.....	39
Cleaning	43
HF cables	45
Cable connections	45
Accessories from other manufacturers.....	47
Sheaths with continuous irrigation and suction	48
Working element 26050 V	51
Optical forceps and scissors	53
Overview of preparation.....	54

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Transvaginal endoscopy, ferti-loscopy, transvaginal falloposcopy 62

GORDTS AND CAMPO set for transvaginal endoscopy	62
Puncture needle 26182 TA	63
Trocar sheath 26182 TC	66
Telescope 26120 BA	68
Sheaths 26182 D/TG	68
Ductoscopes	69
Overview of preparation	71
Examination sheaths, see chapter on 'Hysteroscopes for examination and surgery'	
Flexible forceps and scissors, see chapter on 'Hysteroscopes for examination and surgery'	
Miniature endoscopes, see separate 'Miniature Endoscopes' manual	

Fetoscopy.....72

Working insert with deflection lever	72
Trocars	73
General information	73
Rigid trocars with silicone lip valve	73
Forceps for bipolar coagulation	76
Take-apart forceps	79
CVS biopsy forceps	81
CVS cannula	81
ALKEN micromanipulator	82
Insertion adaptor 26022 Y	83
Overview of preparation	84
Fetoscopy sheaths, see chapter on 'Hysteroscopes for examination and surgery'	

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**Amnioscopes, vaginoscopes,
cystoscopes, urethrosopes,
proctoscopes85**

Micro blood examination set.....	85
Amnioscopes	88
Proximal light deflectors.....	88
Forceps with spring handle, sponge forceps, foreign body/IUD grasping forceps	89
Knife holder 26026	90
Catheter-guiding forceps	91
Universal cystoscope.....	92
Illumination rod	93
Overview of preparation.....	95
Cystoscope/urethroscope sheath, <i>see chapter</i> <i>on 'Hysteroscopes for examination and surgery'</i>	
Semirigid instruments, <i>see chapter on</i> <i>'Hysteroscopes for examination and surgery'</i>	
Biopsy forceps 26958 Z, <i>see chapter on</i> <i>'Hysteroscopes for examination and surgery'</i>	

Maintenance96

Operating/storage conditions	96
Disposal	96
Directive compliance	96
Limitation of liability.....	96
Warranty.....	97
Important information	97
Servicing and repair.....	98
Repair and replacement program	98

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Using this manual

This manual does not contain instructions on how to use the instruments, but instead provides general information regarding assembly/disassembly, cleaning/disinfection, care and sterilization of gynecology instrument sets from KARL STORZ.

This manual presumes knowledge of the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' manual, in which the various cleaning and sterilization methods are described in detail. If you do not yet have this manual, it can be requested without charge from KARL STORZ.

The section 'General information on cleaning and disinfection' contains fundamental instructions, knowledge of which is essential for correct cleaning, disinfection and sterilization.

Instruments whose cleaning, dismantling and reassembly as well as disinfection and sterilization do not require any special knowledge are not listed separately in these instructions.

Particularly sensitive product groups, such as telescopes or fiberscopes, have their own manuals. Separate instruction manuals and service manuals are available for units. Other general regulations (Normal use, User qualification, Safety rules when using instruments, etc.) must be heeded.

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WARNING: Risk of injury/material damage if the instruction manual/warnings are not observed: Read and follow the warnings and instruction manuals for the products and the products used in combination.

WARNING: Risk of injury to the patient/operator and/or damage to the products: Combinations of medical devices are only assured to be safe if they are identified as such in the respective instruction manuals or the intended purpose and interface specifications of the devices used in combination permit this.

WARNING: Risk of injury for patient/user: Damaged or incomplete instruments/accessories may not function correctly. Always check that instruments/accessories are intact and undamaged before and after every use. Exchange incomplete/damaged instruments/accessories and do not use further.

WARNING: Risk of injury for patient: The incorrect fitting of an instrument's assembled components can result in injury to the patient. Check the secure fitting of the components before each application by pulling them gently and/or turning them in both directions without exerting excessive force.

WARNING: Risk of infection for patient: These instruments are not sterile when delivered. Clean, disinfect and sterilize the instruments before initial use and before every subsequent use using validated procedures.

WARNING: Risk of infection for patient/user: Inspect the instrument for visible soiling. This indicates that preparation has not been performed or has been performed incorrectly. In this case, preparation and sterilization must be repeated using validated procedures.

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WARNING: Risk of infection by improper prepa-
ration: The 'Cleaning, Disinfection, Care and
Sterilization of KARL STORZ Instruments' manual
must be followed.

WARNING: Risk of infection and damage to
the medical devices: Effective preparation can-
not be guaranteed if non-validated procedures
are employed. The list of materials approved for
preparation must also be observed. This can be
found on the Internet (www.karlstorz.com) and in
the 'Cleaning, Disinfection, Care and Sterilization
of KARL STORZ Instruments' manual.

WARNING: Risk of infection for patient/user: The
successful sterilization of a medical device can
only be ensured when preparation is performed
under aseptic conditions.

WARNING: The listed cleaning, disinfection and
sterilization methods alone do not guarantee steri-
lity. This can only be achieved if staff follow
recognized and validated preparation methods.

WARNING: For all cleaning of contaminated
endoscopes and instruments, the Accident
Prevention Regulations (UVV) (e.g., wearing of
protective gloves and safety goggles, etc.) must
be observed.

CAUTION: The following information on cleaning
and sterilization only applies to the thermostable
instrument set described in this manual. Special
rules apply for telescopes and fiberscopes (see
manuals for telescopes and fiberscopes).

CAUTION: When preparing and using the solu-
tions, follow the chemical manufacturers' instruc-
tions paying close attention to proper concentra-
tion and exposure times. A maximum immersion
time of 60 minutes (including in water) must not
be exceeded. Prolonged immersion can lead to
damage. After cleaning it is recommended to rinse
the instruments with water.

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CAUTION: Constant changing between various preparation methods is especially hard on the materials and should therefore be avoided. However, a one-off change of method, resulting from the procurement of new cleaning and sterilization equipment, for example, is possible.

WARNING: Risk of injury for patient/user: The instrument may only be used in accordance with the medical rules and procedures for endoscopy.

WARNING: Risk of injury for patient/user: Overloading the instrument or accessory (as a result of applying or exerting too much force on it) may cause the medical device to break, bend and malfunction and injure the patient or user as a consequence. Never overload the instrument. Do not bend bent instruments (xxx, xxx, etc.) back to their original position.

WARNING: Risk of injury for patient! The instruments may only be used by persons with an appropriate medical qualification and who are acquainted with the technique.

WARNING: Risk of burns for patient/user! The high level of light energy means that the distal end, the light port and adjacent components may become so hot that they can cause burns to the patient/user and/or surgical accessories! Failing to maintain sufficient distance between the light outlet and tissue can result in tissue damage.

Avoid direct tissue contact with the distal end as well as with the light port of the telescope and the light cable.

Never rest the endoscope or the light cable on the patient or where it is in direct contact with operating accessories.

Always select the lowest light output of the cold light source which still allows optimal illumination of the operating field.

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WARNING: Danger of glare! Never look directly into the open end of a connected light cable or endoscope.

CAUTION: The sealing systems used in combination with these instruments are recommended for single use only. The efficacy of the sealing systems is no longer guaranteed if they are used more than once.

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**General instructions on preparation
of medical devices**

The medical device cycle shown below illustrates
the sequence of steps for advisable instrument
preparation:

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According to the instrument cycle, preparation of the instruments described here generally comprises the following 6 steps:

CAUTION: The instructions for cleaning and sterilization in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' must be observed. The methods for cleaning, disinfection and sterilization are explained therein in detail.

NOTE: Recommended cleaning accessories and sterilization containers can be found in the catalog chapter 'Cleaning, disinfection, care and sterilization'.

1. Manual preliminary cleaning

The instruments may require manual preliminary cleaning on account of their design (e.g., slits and long, narrow lumens) and depending on the type and degree of soiling.

In addition, a prolonged waiting time before preparation (which allows the contamination to dry on) may make preliminary cleaning necessary while still in the operating room. Organic residues and pharmaceutical products, for example, must be removed mechanically (by brushing, rinsing through with a water gun, etc.) and the instruments must subsequently be rinsed off with completely demineralized or microbiologically pure/sterile water.

Used and contaminated instruments should always be prepared with the minimum possible delay.

Separable instruments must be disassembled following the instructions. Articulating instruments and scissors must be treated whilst open.

Tissue remnants sticking or baked to the jaws must be removed with a cleaning brush, and the instrument must be pre-cleaned with a 3% H₂O₂ solution as soon as possible after use. Use a disinfected cleaning brush to remove tissue remnants sticking or baked to the jaws. Special care must be taken when cleaning the jaws, as they are precision-engineered parts.

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CAUTION: When preparing and using the solutions, the manufacturer's instructions concerning dilution and immersion time must be strictly adhered to. Prolonged immersion can lead to damage.

NOTE: A list of products approved by KARL STORZ for manual cleaning/ disinfection can be found in the appendix of this manual and on the Internet (www.karlstorz.com).

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Thoroughly clean surfaces, interiors and jaws with a suitable soft disinfected brush, sponge or cleaning gun, and rinse. Take care when cleaning small parts.

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After cleaning, rinse the instruments copiously with microbiologically pure/sterile water and dry with sterile compressed air.

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As there is sometimes strong movement of the water, it is necessary to secure all parts in holding devices to avoid damage. To ensure that the instrument is adequately rinsed out, it is essential to connect up the hollow spaces appropriately to irrigation tubes or irrigation devices.

It is not advisable to use lubricants etc. in the machine because they may cause deposits to form on the instruments and may damage plastics. Selective manual care and lubrication of joints etc. is safer and more effective.

Adequate final rinsing with completely demineralized water is required to remove all chemical residues.

① NOTE: A list of products approved by KARL STORZ for machine cleaning/ disinfection can be found in the appendix of this manual and on the internet (www.karlstorz.com).

4. Functional testing

All instruments should be checked for correct functioning before sterilization, at the latest. This functional check should include the following aspects in particular:

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Damaged instruments must be separated out and sent for expert, qualified repair. For this purpose we recommend the KARL STORZ Repair Service, which ensures consistent quality and guaranteed safety of your instruments and units.

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① NOTE: Apply a thin coat of special grease to sealing and running surfaces, such as stopcocks and plungers (Art. no. 27657). This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the stopcocks leakproof. Oil is not suitable as a sealant.

⚠ WARNING: Only use special silicone- or oil-based lubricants, as the sterilants cannot penetrate all silicones and oils.

5. Sterilization

For almost all instruments, steam sterilization by the fractionated vacuum method at 134°C, gas sterilization (formaldehyde, ethylene oxide), plasma sterilization and chemical sterilization (STERIS) are possible. Steam sterilization at 134°C is usually the gentlest and most economical method of sterilization. This also applies to rigid telescopes because the longer exposure time for steam sterilization at 121°C can have a negative effect.

Sterilization is only ever possible on surfaces that have been cleaned and disinfected.

Almost all instruments can be sterilized without being dismantled. Only tighten screws and closures gently, so that they are easy to undo later. Insert the valves into the stopcocks and tighten the screws only gently before sterilization.

① NOTE: If the design of the instrument necessitates sterilization in the dismantled state to ensure effectivity, this is expressly stated for the article in question in this manual.

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For steam sterilization (with saturated steam at 134°C + 3°C), the process defined in EN 285 or DIN 58946 Part 1, Number 3.25.2 should be used. Use prevacuum methods for instruments with a lumen. The steam must fulfill the requirements of EN 285, Appendix B, in order to avoid discoloration and stains.

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CAUTION: When sterilizing telescopes, please follow the instructions in the Telescopes manual.

WARNING: Long lumens, slits and channels, and assembled instruments can only be sterilized completely with the fractionated prevacuum method.

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Refer to the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ instruments' for more detailed information.

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Suitable accessories, storage systems, and sterilization trays for special instrument sets can be found in the relevant chapter of the TELEPRESENCE catalog.

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
Hysteroscopes for examination and operation

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- When cleaning sheaths pay special attention to the sheath lumen and stopcocks.
- For sterilization following cleaning, the stopcocks must be lubricated and disassembled.

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Clean the shaft lumen with a flexible brush.

If the instrument has a working channel: remove sealing caps, and clean with a disinfected flexible brush. Rinse with water and dry with sterile compressed air.

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Clean the inlet and outlet nozzles with a disinfected brush; if necessary, use a cleaning gun with water and sterile compressed air.

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Disassemble the stopcocks for further cleaning.
Unscrew and remove the spring cap, and take
out the cock plug.

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Clean the cock plug thoroughly with a brush.

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
Clean the inlet and outlet nozzles with a brush.

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Rinse out with water and dry with sterile
compressed air.

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Assembly

Lightly grease the contact surfaces on the cock plug ② (special grease for stopcocks; Art. no. 27657).

Note: This thin layer of grease not only serves to lubricate the metal-to-metal surfaces in contact with one another, but also prevents the stopcocks from leaking. Oil is not suitable as a sealant.

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Insert the cock plug and fix with the spring cap. Tighten slightly.

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Take the cock plug apart again for sterilization.

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Semirigid instruments

Thoroughly clean the jaws with a brush.

Then clean the instrument in an ultrasound bath


Rinse the instrument with water. Then dry with sterile compressed air.

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
Place a drop of instrument oil (Art. no. 27656 B) in the joints, holding the instrument in such a way that the oil can flow in and also wet the wire pull.

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
Flexible instruments

Clean the jaws thoroughly with a disinfected brush.

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

Undo the thumbscrew ① and remove the rear handle part from the groove of the forceps insert.

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
Pull off the front handle part toward the rear.
Then clean in an ultrasound bath. Rinse out with
water and dry with sterile compressed air.

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Then lubricate the joints with a drop of instrument
oil (Art. no. 27656 B), holding the instrument so
that the oil can flow into the jacket and thus also
wet the wire pull.

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
Then oil the wire pull from the proximal end.
Afterwards reassemble the instrument in reverse
order.

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
BETTOCCHI polypectomy snare
Disassembly
Unscrew and remove the thumb grip ①.

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
Unscrew and remove the contact socket ②.

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
Hold the handle piece tight, keeping the mount
for the contact socket ② stationary; unscrew the
silver-colored union nut ④ of the guide sleeve
clockwise and pull off toward distal.

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Detach the union nut ④ from the snare insert
by turning counterclockwise, and pull off toward
proximal.

The snare wire, including the plastic covering ⑤,
and lock cannot be dismantled any further, and
must be replaced as a complete unit.

WARNING: Risk of infection! The loop
26159 LC is for single use only and must
not be reprepared.

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Cleaning

Clean the surfaces with a soft disinfected brush. Brush out the hollow spaces. Then rinse the instrument with water and dry with sterile compressed air.

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Assembly

Introduce the new snare insert ⑤ into the guide sleeve and lock by turning counterclockwise.

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Insert the guide sleeve with snare insert into the handle and screw on by turning counterclockwise. When doing this, hold the handle so that the mount for the contact socket ② is stationary.

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Pull the mount for the contact socket ② to proximal. Screw the thumb grip onto the handle. Ensure that the lugs on the handle engage in the corresponding recesses in the thumb grip ①.

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
Screw on the HF contact socket.
Test for proper functioning.

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
Compact miniature hysteroscopes Preparation for cleaning

Remove all instruments from the working channel. Unscrew and remove the light cable and light cable adapter. Unscrew and remove the spring caps (4) and remove the cock plugs. Wipe off grease from the cock plugs. Fully open the (optional) fine metering valve. Brush the lumen of the cock plugs (brush Art. no.: 27650 G). Put the light cable adaptor, cock plugs and spring caps in the small parts tray of the tray holder.

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Hold the instrument inlet by the cylindrical sleeve, and pull off.

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
Pull sealing system(s) off the instrument inlet and insert into the cleaning adapter in the tray holder. Brush the lumen of the instrument inlet (brush Art. no.: 27651 UA). Hold the instrument inlet by the flat part and place onto the side connection of the cleaning adapter.

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
Brush the lumen of the irrigation channel inlet and outlet.

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
Brush the lumen of the working channel ⑤ (brush Art. no.: 27651 UA).

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Assembly

Hold the instrument inlet by the flat part and slide it onto the mount ⑤.

① **NOTE:** If the instrument inlet does not click into place (small picture), slide the cylindrical sleeve to the flat side and then assemble.

Insert both cock plugs ④ into the irrigation ports ③+⑥ and fix with spring caps ② (tighten gently).

Screw the light cable adapter onto the light inlet piece ①.

Test the functioning of the instrument.

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
Preparation for cleaning

1. If necessary, remove instruments from the working channel.
2. Unscrew the light cable and light cable adaptor.
3. Remove the valves for cleaning (when not in the locked position) from the handle by pulling sharply.
4. Inspect all O-Rings for damage and replace if necessary.

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5. Replace damaged sealing rings.
6. Remove the outer sheath ① by holding the working insert ② on the union nut and pulling it forward.
The outer sheath ① engages with a clicking sound and can be pulled off toward the front.
7. Inspect the O-Ring from the working channel and replace if necessary.
8. Remove the sealing cap ⑦ from the compact hysteroscope ②, check again and replace if necessary.

The following **replacement seals** can be found in the spare parts box, **art no. 26252SP**:

10 x 5921100 O-Ring (for valve 26252BV)
10 x 5922700 O-Ring (for valve 26252BV)
10 x 5936400 O-Ring (for valve 26252BV)
5 x 7919291 Sealing cap (for working channel)
5 x 5922000 O-Ring (for outer sheath 26252BO)

9. Press the locking mechanism on the handle and remove the Monobloc adaptor ④ from the valve connector ③.
10. Check the O-Ring of the Monobloc adaptor ④, and replace if necessary.
11. Place the light cable adaptor, the sealing cap (7919291) and Monobloc adaptor ④ in the small parts tray of the tray holder for the final machine cleaning and disinfection.

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Manual cleaning and disinfection Outer sheath (26252BO)

1. Brush through the outer sheath (26252BO) thoroughly and rinse it through with microbiologically pure/sterile water.
2. Then dry with sterile compressed air.

Compact hysteroscope (26252BB/BH)

1. Thoroughly clean (with a disinfected brush) on the handle of the B.I.O.H. shaft all (lumen)
 - both lumina of the Monobloc connector ④ (brush art. no.: 27650 B),
 - the lumen of the working channel ⑥ (brush art. no.: 27651 UA),
 - both lumina of the B.I.O.H. handle which are visible when the suction and irrigation valves ⑧ are removed.

NOTE: The cleaning adaptor 39501XCA can also be used for the manual cleaning of the B.I.O.H. compact hysteroscope.

2. Place the cleaning adaptor 39501 XCA in the suction and irrigation ports (with suction and irrigation valves removed) and thoroughly rinse the lumen of the suction and irrigation connectors using a rinsing aid (e.g., syringe).
3. Rinse through (with microbiologically pure/sterile water) on the handle of the B.I.O.H. shaft all (lumen)
 - both lumina of the Monobloc connector ④
 - the lumen of the working channel ⑥
 - both lumina of the B.I.O.H. handle which are visible when the suction and irrigation valves ⑧ are removed.
4. Then dry all parts (with sterile compressed air) and on the handle of the B.I.O.H. sheath all (lumen)
 - both lumina of the Monobloc connector ④
 - the lumen of the working channel ⑥
 - both lumina of the B.I.O.H. handle which are visible when the suction and irrigation valves ⑧ are removed.

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
Care

1. Oil the suction and irrigation valves ⑧: apply a drop of instrument oil (art. no. 27656 B) to the opening (bistable valve) or between the top and the metal mount (monostable valve).

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
1. Insert the cleaning adaptor 39501 XCA into the valve holes.
2. Secure the valves in an open position in the tray holders.
3. Position the light cable adaptor on the specially designed holders.
4. Position the tubes of the tray holder on the cleaned Monobloc 26252 BM.
5. Place the BETTOCCHI compact hysteroscope correctly orientated in the tray holder.
6. To do this, open the outer sheath. This, however, remains on the compact hysteroscope according to its fixed position in the tray holder.
7. Place the light cable in the holders on the wall of the tray holder, put the cover onto the tray holder and lock into place.
8. Place the equipped tray holder into the washer and disinfectant and connect the LUER connectors on the front side of the tray holder.
9. Prepare the instrument set in the washer and disinfectant with thermal disinfection at up to 93°C (device-specific).

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
**Hysteroscopes for
examination and
operation**

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- ① Telescope (without fitted sheath) 26008 BAC
- ② Telescope 26008 BAC with fitted surgical sheath ④ 26152 DB (position 2).
- ③ Telescope 26008 BAC with fitted examination sheath ⑤ 26152 DA (position 1)
- ⑥ Sealing cap 6809295 for sheath
- ⑦ Sealing cap 27550A for working channel

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Disassembly

Decontamination (cleaning and disinfection) must only be performed when the instrument is disassembled.

To this end, remove the telescope by pressing the locking mechanism on the sheath and pulling it out of the sheath.

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Manual cleaning and disinfection

CAUTION: Telescopes must not be cleaned in an ultrasound bath under any circumstances.

- Clean the surfaces carefully with a disinfected sponge or a disinfected, soft brush. Brush the lumen through thoroughly.
- Only cleaned and disinfected instruments can be sterilized effectively.
- Immerse all the parts in a disinfectant solution making sure that all air bubbles have been removed and that the parts are completely covered, and taking the chemical manufacturers' instructions on concentration and exposure times into account.
- Rinse thoroughly with completely demineralized water and then dry with sterile compressed air.

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Manual and machine preparation (Tristel)

NOTE: Tristel: Nationally applicable regulations concerning the use of Tristel must be observed.

- Manual/machine preparation using the Tristel method is permitted (up to 200 preparation cycles)
- Exceeding the 200 cycles may lead to material changes. Check the instruments for damage before / after every use.

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Machine Cleaning/Disinfection

Ensure that all parts are stored correctly to avoid any damage occurring when machine cleaning is in progress. To this end, the parts are best placed in a tray holder.

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Sterilization

Cold sterilization with Tristel (up to 200 cycles) may only be performed if the respective national/ official regulations are observed. Familiarize yourself with these regulations and ensure that they are strictly observed.

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Care

Oil moving parts locally at the contact points with instrument oil (Art. no. 27656). Test for proper operation.

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Assembly

Insert the telescope into the sheath By pressing the locking button, the sheath can be moved on the telescope at will.

To lock in position 1 or position 2, release the locking button. The sheath is now locked in the required position.

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References:

HF loops and electrodes

see chapter on 'Resectoscopes in intrauterine HF electrosurgery'

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Safety instructions for the use of resectoscopes

The working element, sheath, electrodes and high frequency cable must be examined for damage each time they are used. In particular, the high frequency cable must be carefully examined for breaks or faulty plug-in connections.

CAUTION: Sheaths with ceramic tips must be handled carefully. The sheaths must be checked for cracks or other damage each time they are used. Damaged sheaths (chips or cracks) must never be used, since malfunctions could occur which could lead to patient injury. Only use obturators specially designated for sheaths with ceramic tips. Unsuitable obturators may cause damage to the ceramic tip.

WARNING: Damaged HF cables and plug-in connections may result in serious danger for the patient and/or user and must not be used under any circumstances.

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
**Model 26040 xx, 26050 xx and
26055xx**
(For use with electrodes with 2 guide rods)

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To change the electrodes, press the side locking button ① (locking function) and pull out the electrodes toward the front.

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Clean the working element with a disinfected sponge or brush.

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
Flush through with water and dry with sterile
compressed air.

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
The bare, distal electrode parts must not be bent
or deformed. The correct configuration of the free
electrode end is essential for correct function-
ing and safety in regard to the electric strength
defined for the system.

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Damaged electrodes or ones which have become
thin because they have burned down must be
replaced.

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WARNING: It is essential to read the
safety information in the chapter on HF
electrodes as well as the instruction
manual for your HF unit before using the
resectoscope.

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
Model 26040 xx and 26055 xx
(For use with electrodes with 1 guide rod)

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To change the electrodes, press the side locking
button ① (locking function) and pull out the elec-
trodes toward the front.

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
Clean the working element with a disinfected
sponge or brush.
Clean the lumen with a flexible brush.

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
The bare, distal electrode parts must not be bent or deformed. The correct configuration of the free electrode end is essential for correct functioning and safety in regard to the electric strength defined for the system.

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Damaged electrodes or ones which have become thin because they have burned down must be replaced.

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WARNING: It is essential to read the safety information in the chapter on HF electrodes as well as the instruction manual for your HF unit before using the resectoscope.

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HF electrodes: Safety instructions

Incidents in high frequency surgery which cause damage to the instrument set or even injuries to the user or patient are virtually always attributable, on closer investigation, to faulty application or defective parts of equipment. The chief sources of error prove to be damaged cables and/or electrodes. Another cause may be the use of unsuitable different makes of accessories which are not compatible with KARL STORZ equipment.

To prevent damage or hazards, it is imperative that users acquaint themselves with the following safety information.

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Configuration of the loop/hook electrode

Damaged insulation on HF electrodes may become a source of danger. The high peak voltages delivered by commercially available HF-generators impose stringent demands on the electrical insulation of electrodes. Inspect electrode insulation for evidence of damage or current leakage prior to use and discard any suspect electrodes. With resectoscopes, it is particularly the correct distance between the loop and the endoscopic telescope which is important. Safe and perfect operation of a resectoscope essentially depends on the loop being located at a precisely defined distance from the end of the sheath and from the conductive metal tip of the endoscopic telescope. The structurally determined clearance and creepage distances determine the maximum permissible peak voltage (U peak) for the entire system. Any change in this distance can lead to a voltage flashover and thus to a flow of current from the electrode to the endoscopic telescope and then to the sheath. Any current flowing via the shaft can in turn result in burns to nearby tissue, in particular if other parameters, for example lubricants, have an influence on the current path.

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Resectoscopes for intra- uterine HF Electrosurgery

Use of lubricants on metallic sheaths can readily lead to electrical hazards, since, in spite of their low electrical conductivities, lubricants can locally break down and transmit high current densities. Damage to the instrument set, particularly telescopes, can also be caused by sparks at the tip of the endoscope.

From experience, bent loops are no seldom occurrence in everyday clinical situations. Knocking the loops against hard objects (e.g. instrument table) or careless removal of adhering resected material from the loops may tend to bend the loops. Loops require careful handling and storage. KARL STORZ therefore offers special protective sleeves for storing electrodes. Intentional bending of the loops by the user to achieve a supposedly improved incision or to improve the quality of incision with older loops is just as hazardous. Such actions mean that the loop no longer has the configuration necessary to ensure that the instrument is handled safely. Loops where the wire has burned down considerably in comparison with new loops must be discarded. Such loops may break if you continue using them. The patient and user can be placed at risk from this, and other instruments may be damaged.

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Application

WARNING: The instructions and interface specifications for medical devices used in combination must be observed precisely.

WARNING: Electrodes labelled as unsterilized must be thoroughly cleaned, disinfected and/or sterilized before being used for the first time and before each subsequent use.

WARNING: Ensure that the electrodes are in perfect condition before every use. The following points must be checked carefully in particular:

- Is the insulation undamaged?
 - Are the HF electrodes and contact ends undamaged and in perfect condition?
- If there is any sign of damage, the electrode must not be reused under any circumstances.

WARNING: As a rule, HF electrodes in combination with standard resectoscopes must only be used with a recovery peak voltage of 2.0 KVP max. in the standard cutting and coagulation modes. System combinations which permit lower or higher peak voltages are separately identified.

WARNING: The electrodes must be fully inserted into the working element until they click into place. Electrodes which have not been inserted properly can cause disruptive discharges into the working element.

WARNING: Only activate the HF current when the electrode is in your field of vision and the required tissue contact has been made.

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Resectoscopes for intra- uterine HF Electrosurgery

WARNING: Do not activate the HF current if the HF electrodes are in an air or gas bubble (e.g. at the roof of the bladder). Danger of explosion!

WARNING: Only use the power setting absolutely necessary for the intervention. If the required effect is not achieved in spite of a standard setting on the HF unit, all components involved must be carefully checked for possible damage and correct connection before increasing the power output.

WARNING: When using electrosurgical equipment, other equipment may suffer electromagnetic disturbance, particularly cardiac pacemakers. For this reason, preventive measures should be taken to ensure the patient's well-being in case of such a disturbance.

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Cleaning

HF electrodes should be wiped off during the operation between uses with a moist compress to clean the coagulating surfaces. Incrustations can be dissolved by briefly immersing in a 3% H₂O₂ solution. Then rinse off with sterile water and wipe off again.

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Manual cleaning

Carefully remove tissue remnants from the electrodes with a brush.

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Sterilization

! **CAUTION:** HF electrodes must not be disinfected in solution. Moisture or residue of disinfectant/cleaning agents on the contact connections may lead to malfunctions during operation.

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
The electrodes must be placed in a protective sheath for sterilization and storage (Art. no. 280 and 280 L).

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Inserting the electrode

Slide the electrode stabilizer tube over the tube of the working element. Push the contact end of the electrode into the opening provided on the working element until it clicks into place.

NOTE: Color markings on electrodes and sheaths make identification easier. Only use electrodes and sheaths with the same color marking. Parts coded with two colors are compatible with both.

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Manual cleaning and disinfection

HF cables can be cleaned with a cleaning agent/disinfectant specially recommended for endoscopes and approved by KARL STORZ, using a disinfected brush or sponge.

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Preparation by machine

HF cables can be prepared by machine both chemically and thermally using special programs for endoscopic instrument sets.

NOTE: We recommend using water at all stages of the program. The selection of chemicals and their concentration must also be observed precisely. Lubricants should not be used in the machine because coatings may form on the metal contacts and plastics may be damaged.

CAUTION: After immersing the cable in solution carefully dry it out in order to avoid any malfunctions due to residual moisture at the mating contacts. After cleaning, steam sterilization should preferably be performed at 134 °C using the fractionated prevacuum method.

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WARNING: The HF cable is delivered unsterilized and must therefore be cleaned and sterilized before the first as well as each subsequent use.

NOTE: Make certain that the plug connections maintain good contact. Pay attention to locking systems or markings on the plugs.

WARNING: Poor or insufficient contact may lead to voltage flashovers due to the high voltage at the contacts, which could damage the contacts or the insulation in the plug.

WARNING: HF cables are subject to a certain degree of wear. Damage to the cable insulation or to the plug connections may lead to HF voltage flashovers to conductive metal parts. Parts thus unintentionally at HF potential, such as an instrument, could become dangerous for the patient and user and possibly lead to burns. It is therefore essential to ensure that the HF cable and connectors are in a perfect condition each time they are used. Special attention must be paid to the following points:

- Burn marks or traces of flashover on instrument connectors
- Defects on the plug connections
- Breaks in the cable or cable insulation

WARNING: If there is any sign of damage, the HF cable must not be reused under any circumstances.

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**Accessories from other
manufacturers**

Safe use of KARL STORZ instruments is only assured with KARL STORZ high frequency cables and electrodes. The use of other makes of cable which are not compatible with KARL STORZ instruments may lead to the hazards mentioned in the section about cable connections.

Only the following KARL STORZ HF cables may be used for KARL STORZ resectoscopes:

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
Inspections have indicated that a not insignificant number of cutting loops which were supplied by competitors and explicitly described as compatible with KARL STORZ resection sheaths were in fact incompatible. Too small a distance between the loop and the tip of the endoscope, as described in the section on electrodes, proved to be a major source of danger to the patient and user. KARL STORZ can accept no liability for any damage or injury resulting from the use of unsuitable accessories supplied by competitors or damaged accessories.

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**Sheaths with continuous irrigation
and suction**
Inner sheaths, fixed or rotatable
The sheath consists of an inner and outer
sheath which are cleaned individually.

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Clean the sheath lumen using a suitable disin-
fected flexible brush.

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If necessary, a water and compressed air gun
may be used.

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Clean the nozzle first with the stopcock open,
(cock plug parallel to the nozzle) using a suitable
flexible brush.

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Then rinse out the stopcock. Rinse until clear
water runs out the front end.

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The stopcocks must be disassembled for further
cleaning. Loosen the spring cap on the stopcock,
and remove cock plug.

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Clean the cock plug with a disinfected brush.

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
Also clean the inlet and outlet nozzles with a
brush.

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The nozzle and the sheath lumen must be rinsed
out well on the inside each time cleaning is per-
formed. Then dry using sterile compressed air.
Again, the cleaning gun with various attachments
may be used for this.

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Insert the cock plug and fix in place with the spring cap. Tighten slightly.

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
Working element 26050 VN
for use with GYNECARE Versapoint
snare wires and electrodes

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- ① Seal
- ② Guide groove
- ③ Sheath lock
- ④ Mount for HF electrode
- ⑤ Telescope lock

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Disassembly for cleaning
Pull the electrode on the cable connector side-
ways out of the mount ④ and remove.

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Pull off the seal ①.

① **NOTE:** The seal is a disposable item and
must be replaced with a new one after
each use. Clean the working element with
as soft a disinfected brush as possible.

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Clean the lumen with a disinfected flexible brush.

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To clean the snare mount, remove the heavy dirt with a soft brush, then immerse it in a cleaning and disinfecting solution. If necessary, clean again with a brush or in a machine.
Oil the moving parts.

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Assembly

Install a new seal in the correct position ①.
Insert the electrode in the lateral guide groove ② and press it into the mount as far as the limit stop ④.

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Optical forceps and scissors

Clean the exterior of the forceps with a sponge,
and the jaws with a disinfected brush.
Brush out the telescope channel.
Next, clean in an ultrasound bath
Rinse out with water and blow dry with sterile
compressed air.

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Oil the joints of the forceps jaws, holding the
instrument in such a way that the oil is able to
flow from the distal end into the joint and wet the
wire pull.

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Press and hold the locking button ⑤ (do not turn).
The handle is now opened completely in order to push the sheath tube ② out towards the front.


CAUTION: Do not turn locking button ⑤, as this may cause it to snap off.

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To remove the internals, hold the jaws ① tightly and open the bayonet mount ⑦ by rotating the sheath tube ② ¼ turn. Pull out the internals.

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CAUTION: Please note that in CLICK'line instruments in the 2 and 3 mm sizes and in curved instruments the internals cannot be removed from the sheath tube.

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
Cleaning

NOTE: For heavy soiling, we recommend using an ultrasound bath for a maximum of five minutes.

Manual cleaning

Remove tissue residue with a cleaning brush. Brush and rinse out sheath tubes on the inside. Blow dry with medically clean compressed air.

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

Machine cleaning

Make certain all parts are positioned securely so that damage during cleaning is prevented. The sheath tube must be rinsed out on the inside. By means of the LUER-lock irrigation connection ④, easy attachment to a syringe or cleaning machine is also possible.

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
Assembly

Check all parts for damage. Replace if necessary.
Check jaws and joints especially carefully, and also the insulation on instruments for HF use. Insert the internals from the front into the sheath tube ②.
Lock the internals by rotating ¼ turn in the bayonet mount on the sheath tube ②.

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
Place a drop of instrument oil (Art. no. 27656 B) in the joints.

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
In the case of handle 33140, the thumb grip can be removed. For this purpose, pull the thumb grip toward distal and remove by screwing in clockwise direction.

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
Uterine forceps

Attached cones can be pulled off toward distal to clean the instrument. Then clean all outer surfaces with a disinfected sponge.

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
Thoroughly clean jaws and cones with a disinfected brush.

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Clean the perturbation channel with a brush and rinse out with a cleaning gun. Then clean the instrument in an ultrasonic bath, rinse out with water and dry with sterile compressed air.

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Oil the joint (instrument oil Art. no 27656 B) while holding the instrument in such a way that oil can flow into it.

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
**COHEN intrauterine probe
26168 UN**

The cones can be pulled off toward distal to clean the instrument.

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Brush out the lumen of the sheath and cones. Then clean the instrument in an ultrasonic bath, rinse out with water and dry with sterile compressed air.

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

To remove the adjusting sleeve, turn the thumb-screw clockwise until the sleeve can be pulled off toward distal.

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Cleaning

Dismantle and clean the stopcock (see Sheaths with Luer lock stopcock).

Brush out the lumen of the sheath and cones. Clean in an ultrasound bath. Rinse and dry with sterile compressed air.

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Electrode handle 265200 43

The electrode handle can be cleaned in the washing machine at up to +95°C or in an ultrasound bath.

Then rinse the connector and loop mount with water and dry with sterile compressed air.

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**GORDTS and CAMPO set for
transvaginal endoscopy**

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- ① Release switch
- ② Tensioning lever
- ③ Union nut
- ④ Limit stop
- ⑤ Scale for limit stop
- ⑥ Scale for dilation sheath
- ⑦ Knurled nut
- ⑧ Limit stop for trocar sheath
- ⑨ Lock, open
- ⑩ Thumbscrew
- ⑪ Cock plug
- ⑫ Spring cap
- ⑬ Light connection piece
- ⑭ Telescope lock

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
Puncture needle 26182 TA

Assembly

Before commencing assembly, actuate the release switch ① to ensure that the spring is not under tension.

WARNING: If the spring is under tension during assembly, there is a danger that the needle will suddenly jump out and cause injury!

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
Slide the dilation sheath 26182 TB over the puncture needle 26182 TA and lock in place.

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
Slide the trocar sheath 26182 TC over the dilation sheath until it reaches the mechanical stop. Make sure that the lock is open during insertion.

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
Tighten the tensioning lever (2) until you feel and
hear it click into place.


WARNING: As soon as the needle
has been loaded, it can be unlocked by
pressing the release switch.
Danger of injury!

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
Set the required penetration depth of the punc-
ture needle 26182 TA. Adjust the limit stop (4)
to the required depth by pressing both surfaces
simultaneously, and engage.

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
Disassembling the puncture needle 26182 TA

Adjust the limit stop (4) to a penetration depth of
25 mm by pressing the two surfaces of the limit
stop (4) simultaneously.

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Release the puncture needle by pressing the
release switch (1).

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

Unscrew the union nut ③ and remove the inner part with the spring toward proximal.

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
Remove the limit stop ① by simultaneously pressing both surfaces.

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Caution: The needle may only be used once. When reassembled, a new needle from pack 26182 TAA must be used.

To install the needle, proceed in the reverse order.

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Cleaning

Brush out the hollow spaces. Next, clean in an ultrasound bath. Rinse with water and dry with sterile compressed air.

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
Trocar sheath 26182 TC


Disassembly

Remove the sealing cap. Unscrew and remove
the spring cap.

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
Remove the cock plug ⑪.
Thoroughly clean the cock plug and the collar
using a brush.

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Unscrew and remove the trocar head by releas-
ing the thumbscrew ⑩.

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Take out the seal ① and seal holder ②.
Cleaning: Inspect the seal and, if necessary,
replace with a new one.

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
Remove and inspect the silicone lip valve ③; if
necessary, replace with a new one.

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
Clean the sheath lumen with a soft disinfected
brush.

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Rinse with water; dry with sterile compressed air.

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Also brush the trocar head. Rinse with water; dry
with sterile compressed air.

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Assembly

To assemble, proceed in the reverse order. Lightly grease the contact surfaces on the cock plug ⑪ (special grease for stopcocks; Art. no. 27657).

⑪ **NOTE:** This thin layer of grease not only serves to lubricate the metal-to-metal surfaces in contact with one another, but also prevents the stopcocks from leaking. Oil is not suitable as a sealant.

Test the cock plug for proper operation.
Take the cock plug apart again for sterilization.

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Telescope 26120 BA

⑪ **NOTE:** More detailed instructions on the cleaning and care of telescopes can be found in the HOPKINS® Telescopes Manual.

Unscrew both light connection adapters before cleaning.

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Sheaths 26182 D/TG

Remove heavy soiling using a soft brush. Brush out the lumen, then rinse with water and dry with sterile compressed air.

For cleaning the LUER lock stopcocks, see chapter on »Sheaths with LUER lock connector«

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
Ductoscopes Preparation for cleaning

NOTE: More detailed instructions on the cleaning and care of telescopes can be found in the HOPKINS® Telescopes Manual.

Immediately after use, clean the exterior of the ductoscope of heavy soiling (secretions, etc.). Use a disposable cloth soaked in active disinfectant solution. Brush through the channels and rinse thoroughly with cleaning and disinfection solutions.

CAUTION: Due to the small diameter, it is absolutely necessary to clean the channels immediately after use to prevent debris from drying and sticking which can cause blockages.

CAUTION: The working/irrigation channels of the miniature and neuro-endoscopes should be thoroughly rinsed immediately after use with a syringe (filled with active disinfection solution).

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
- Remove all instruments from the working channel.
- Remove adaptors and other detachable parts from the light inlet ports and the working channel.

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**Transvaginal endoscopy,
fertiloscopy, transcervical
fallopscopy, Ductoscopy**

The glass surfaces of the light carrier must be thoroughly dried with alcohol after cleaning/disinfection. Residues of disinfectant and cleaning agents on the glass surfaces may burn on when the fiberoptic light cable is connected, and this inhibits light transmission considerably. Firmly adhering deposits on the glass surfaces can be removed with special cleaning paste (Art. no. 27661).

References:

Examination sheaths

see chapter on

'Hysteroscopes for examination and surgery'

Flexible forceps and scissors

see chapter on

'Hysteroscopes for examination and surgery'

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Working insert with deflection lever

Remove the obturator.
Clean the sheath with a disinfected sponge or brush.

Brush out the telescope channel. Next, clean in an ultrasound bath

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Rinse out with water and dry with sterile compressed air.

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
Thoroughly clean the distal joint of the working insert with a brush.

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Trocars

General information for trocars and accessories

WARNING: For pointed trocar tongues, we recommend the use of protective caps for protection during transportation and to prevent injuries.

Protective caps may be placed on the individual tongue or onto the complete trocar.

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Rigid trocars with silicone lip valve

NOTE: For disassembly, cleaning and reassembly of LUER lock stopcocks, see the 'Instructions for LUER lock stopcocks' (page 69).

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

To remove the silicone lip valve ①, unscrew the knurled ring ②.

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
Take out the silicone lip valve ①, inspect it for damage and, if necessary, replace it with a new one.

NOTE: We recommend changing the silicone lip valve after each use.

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Assembly

Thoroughly clean the lumen and stopcock with a brush. Insert the silicone lip valve ① into the knurled ring and screw on.
Put on the sealing cap.
Insert the trocar tongue.
Put on the trocar protection cap.

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Instructions for LUER lock stop-cocks

Cleaning before disassembly

NOTE: First clean the nozzle and sheath lumen with the stopcock open; use the cleaning gun if necessary.

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
Open the stopcock (cock plug ② parallel with the nozzle ③) and rinse out. Rinse until clear water flows out of the front end.

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
Disassembly

Loosen spring cap on the stopcock, and remove cock plug.

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Thoroughly clean the nozzle with a brush. The nozzle and the sheath lumen must be rinsed out well inside each time they are cleaned, and then dried with sterile compressed air.

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Assembly

Lightly grease the contact surfaces on the cock plug (special grease for stopcocks; Art. no. 27657).


NOTE: This thin layer of grease not only serves to lubricate the metal-to-metal surfaces in contact with one another, but also prevents the stopcocks from leaking. Oil is not suitable as a sealant.

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Insert the cock plug and fix in place using the thumbscrew. Tighten slightly.

Test the cock plug for proper operation.

NOTE: Take stopcock apart again for sterilization.

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
Forceps for bipolar coagulation
Models 11540 xx

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
- ① Handle
- ② Outer sheath
- ③ Working insert
- ④ Jaw
- ⑤ Unlock button
- ⑥ Telescope
- ⑦ Telescope lock
- ⑧ Jaws

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
Disassembly

Open the telescope lock ⑦: From position 0, turn left 45°. Carefully remove the telescope ⑥ (11540 AA).

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
Turn (11540 FGR, image A) or press (11540 HM, Image B) the unlock button ⑤ into the 'Open' position until it clicks.

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
Pull out the working insert ③ (11540 FGR, Image F/11540 HM, Image G).

WARNING: Hold pointed or sharp jaws with a cloth.

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Manual preparation

Remove tissue remnants with a disinfected cleaning brush. Brush and rinse through the lumen before drying it with sterile compressed air, in particular the area around the plug-in connection.

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
**Machine preparation
Washer/Disinfector**

Ensure that all parts are stored correctly to avoid any damage occurring when cleaning is in process.

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
Assembly

Open handle ①.
Turn (11540 FGR, Image A) or press (11540 HM, Image B) the unlock button ⑤ until it clicks into the 'Open' position.
Insert the working insert ③ (11540 FGR, Image E/ 11540 HM, Bild F).
Turn (11540 FGR, Image A) or press the unlock button ⑤ until it clicks into the 'Close' position. The unlock button ⑤ clicks into place audibly (11540 HM, Image B).

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Open the telescope lock ⑦: From position 0, turn left 45°. Carefully insert the telescope ⑥ (11540 AA).
Close the telescope lock ⑦: From left position, turn upwards 45° to 0.

WARNING: Hold pointed or sharp jaws with a cloth.

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
Take-apart® forceps for bipolar coagulation Models 26167xx, 26184xx


- ① Connecting piece
- ② LUER lock irrigation connection
- ③ Unlock button outer sheath
- ④ Unlock button inner sheath
- ⑤ Working insert
- ⑥ Outer sheath
- ⑦ Thumbscrew
- ⑧ Inner sheath

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Disassembly

Hold the Take-apart® forceps by the jaws of the working insert and by the connecting piece ① behind the handle.
Press the unlock button ③ on the connecting piece ① and simultaneously pull out the working insert toward the front.

WARNING: Grasp pointed or sharp jaws with a cloth.

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
By pressing the unlock button (4), detach the inner sheath (8) and remove toward the front.


Unscrew the outer sheath with the thumb-screw (7), and remove toward the front.

Then remove the sealing cap (9) from the handle, check for damage, and replace with a new one if necessary.

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Cleaning

NOTE: For heavy soiling, we recommend using an ultrasonic bath of five minutes maximum duration.

Manual cleaning


Remove tissue remnants with a cleaning brush. Brush out the lumen, rinse and blow dry with sterile compressed air, especially in the area of the plug-in connections.

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By means of the LUER lock irrigation connection (2), easy attachment to a syringe is possible.


Machine cleaning

Make certain all parts are positioned securely so that damage during cleaning is prevented. By means of the LUER lock irrigation connection (2), easy attachment to a cleaning machine is possible.

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Assembly

Attach sealing cap (9) to handle. Then attach outer sheath to the handle and tighten the thumbscrew (7). Slide inner sheath (8) into outer sheath (6). Hold the inner sheath at the front and slide into the holder in the sliding piece up to the limit. Twist the inner sheath until it clicks into place.

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
Insert the working insert ⑤ with the jaws closed through the outer sheath ④ up to the connecting piece ① and put in the correct position. (Unlock button ③ clicks audibly into the place). The working insert ⑤ must no longer turn in the sheath.

CAUTION: The jaws must not touch when closing.

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CVS biopsy forceps

Clean the exterior with a soft disinfected brush; rinse with water.

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
Rinse out the irrigation channel with water, and then dry with sterile compressed air.

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CVS cannula

Remove guidewire. Brush out the lumen, rinse out with water, and dry with sterile compressed air.

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
- ① Connector for fetoscope sheath
- ② Drive wheel
- ③ Adjusting lever
- ④ Unlock button
- ⑤ Instrument channel

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Cleaning and sterilization

The micromanipulator can be dismantled for cleaning: press the unlock button ④ and remove the insert (adjusting lever with drive wheel) from the micromanipulator housing.

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
After external cleaning, the individual parts can be steam sterilized at up to 134°C or gas sterilized.
Grease the insert and put it back into the housing.

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
Insertion adaptor 26022 Y

Dismantle the insertion adaptor for cleaning;
remove the two seals.

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Inspect both seals and, if necessary, replace
them with new ones.

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Clean the insertion adaptor with a disinfected
brush. Insert the seal and reassemble.

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References:

Fetoscopy sheaths
Obtulators
Flexible instruments
Semirigid instruments
see chapter on
'Hysteroscopes for examination and surgery'
Miniature endoscopes
see Miniature Endoscopes manual

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
Amnioscopes, micro blood examination set, cystoscopes


Micro blood examination set

- ① Instrument in non-extended state
- ② Instrument in extended state (detail)
- ③ Guide pin
- ④ Tip of knife
- ⑤ Capillary tube
- ⑥ Connection for battery light source or for adapter to the fiberoptic light cable
- ⑦ Battery light source
- ⑧ Pull rod (contains scalpel holder)

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

Disassembly

- Remove the battery light source. If you have connected a cold light source as an alternative to the battery light source, remove the light cable and the adapter.
- Remove the disposable capillary tube (5).
- Photos 1 to 4: Press the guide pin (3) down and push it down vertically as far as it will go into the left-hand groove. Then guide it - while still keeping it pressed down - to the right, in order to reach the right-hand groove. Finally, push it upward as far as it will go into the right-hand groove so the pull rod (8) extends out of the instrument.
- Remove the disposable scalpel blade.
- Close the scalpel holder (which is a part of the pull rod). Remove the pull rod from the housing by turning it counterclockwise.

Manual preparation

Carefully clean all surfaces with a disinfected sponge or soft brush. We then recommend the use of an ultrasound bath for a maximum of five minutes. Then rinse with water, in order to remove all cleaning agent residue. Immerse all the parts in a disinfectant solution making sure that all the air bubbles have been removed and that the parts are completely covered and taking the manufacturer's instructions on concentration and exposure times into account. Rinse all parts thoroughly with microbiologically pure/sterile water and dry using a lint-free disposable cloth/sterile compressed air.

Machine preparation

Prior to machine preparation, a manual preliminary cleaning is required (possibly in an ultrasound bath). Ensure that all parts are placed correctly to avoid any damage occurring when machine cleaning is in progress. To this end, the parts should ideally be placed in a tray holder. The use of the following processes is permitted:

- chemical disinfection
- thermal cleaning at 93°C

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

Preparation of the battery light source

Wipe-down disinfection is ideal for the cleaning/disinfection of the battery light source.
Steam sterilization is possible at up to 134°C.

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

Assembly

- Set the guide pin in the uppermost position in the right-hand groove. Twist a new pull rod ⑧ clockwise into the housing. Open the scalpel holder (which is a part of the pull rod).
- Place the disposable scalpel blade into the scalpel holder ⑧ (1st photo). Close the scalpel holder. Press the guide pin down and push it down into the right-hand groove.
- Place the capillary tube into the opening provided (2nd photo).
- Close the battery light source ⑦ or the light cable adapter and the light cable at the other end of the instrument (3rd and 4th photos).
- Move the guide pin from the lowest position in the right-hand groove into the left-hand groove. Then test the instrument for proper operation (see next page).
- Test for proper functioning: If you press the guide pin ③ down and push it from the bottom of the left-hand groove to the top, the tip of the knife will protrude from the instrument. To retract the tip of the knife back into the instrument, press the guide pin down again and push it down as far as it will go into the left-hand groove. Note: The guide pin should only be moved in the right-hand groove when the instrument is to be assembled or disassembled (i.e., not during application).

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Amnioscopes

Remove obturator from the outer sheath.

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Clean exterior surfaces with a disinfected sponge.

Clean sheath lumen with a disinfected brush.

Thoroughly rinse with water and then blow dry with sterile compressed air.

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Proximal light deflectors 26201 H and 26201 HL

CAUTION: Never immerse the magnifier 10338 TA of the proximal light deflector in liquid.

CAUTION: Ultrasound cleaning of the magnifier is not permitted.


CAUTION: Sterilization of the magnifier 10338 TA is not permitted. Magnifier 10338 TA is fully autoclavable.

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Clean and disinfect the magnifier 10338 T of the proximal light deflector using a cloth moistened with disinfectant solution.

Clean the surface of the proximal light deflector, not including the magnifier, using a sponge. Press the spring mechanism and clean the hollow spaces with a brush.

Thoroughly rinse the instrument with water, and dry with sterile compressed air.

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
Unscrew and remove both light connection adapters before cleaning.
The glass surfaces of the light carrier must be thoroughly dried with alcohol after cleaning. Residues of disinfectant and cleaning agents on the glass surfaces may burn on when the fiberoptic light cable is connected, and this inhibits light transmission considerably.
Firmly adhering deposits on the glass surfaces can be removed with special cleaning paste (Art. no. 27661).

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Forceps with spring handle, sponge forceps, foreign body/IUD grasping forceps

Press the guidewire down and disengage from the handle piece.

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
Remove the seal. Take off the LUER lock cap.
Take out the internal part of the forceps toward proximal.
Clean with a sponge or brush.

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
Rinse out with water and blow dry with sterile compressed air.

To reassemble the instrument, proceed in reverse order.
Test the forceps for proper operation.

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Knife holder 26206 Disassembly

Remove the knife holder insert from the sheath.

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Cleaning

Clean the exterior surfaces with a disinfected sponge or brush.
Brush out the sheath lumen.

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Then clean the instrument in an ultrasound bath.

Thoroughly rinse the instrument with water, then dry with sterile compressed air.

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
Assembly

Press the knife holder insert together and insert it into the sheath from proximal. Screw the union nut tight.

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Catheter-guiding forceps 26208

see the section 'Semirigid instruments' in the chapter on 'Hysteroscopes for examination and surgery'

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
Universal cystoscope


Take the rubber seal off the instrument channel.

Unscrew and remove the two light port adapters.

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
The glass surfaces of the light carrier must be thoroughly dried with alcohol after cleaning. Residues of disinfectant and cleaning agents on the glass surfaces may burn on when the fiberoptic light cable is connected, and this inhibits light transmission considerably. Firmly adhering deposits on the glass surfaces can be removed with special cleaning paste (Art no. 27661). See p. 10 for cleaning the LUER lock stopcocks.

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
Clean the exterior with a soft disinfected brush. Brush out the lumen, then rinse with water and dry with sterile compressed air.

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Illumination rod

Before cleaning, unscrew and remove both the
adaptors. The illumination rod can be autoclaved
at up to 134°C, gas or plasma sterilized.

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The glass surfaces of the light carrier must be
thoroughly dried with alcohol after cleaning.
Residues of disinfectants and cleaning agents on
the glass surfaces may burn on when the fiber-
optic light cable is connected, and this inhibits
light transmission considerably.
Firmly adhering deposits on the glass surfaces
can be removed with special cleaning paste
(Art. no. 27661).

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CAUTION: Do not clean the illumination
rod in an ultrasound bath.

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**Amnioscopes, micro
blood examination set,
cystoscopes**

References:
Cystoscope/urethroscope sheath
see chapter on
'Hysteroscopes for examination and surgery'
Semirigid instruments
see chapter on
'Hysteroscopes for examination and surgery'

**Semirigid instruments for the universal
cystoscope**
see chapter on
'Hysteroscopes for examination and surgery'

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Disposal

Standard batteries and rechargeable batteries may contain toxic substances, and must be disposed of in an environmentally friendly way.

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Directive compliance

The medical devices mentioned bear the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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Limitation of liability

As the supplier of these medical devices we only regard ourselves as responsible for the safety, reliability and performance of the products if assembly, system expansions, readjustments, modifications, or repairs have been performed by persons authorized by KARL STORZ, and the product has been used in accordance with its operating instructions at all times.

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Warranty

For **two** years as of delivery to the end customer, we agree to replace the goods free of charge if proof can be provided of faulty materials or workmanship. In doing so we cannot accept the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.
Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
D-78503 Tuttlingen/Germany

Any repairs or modifications to medical devices by unauthorized persons shall release us from any liability for their safety. Any such repair or modification performed during the warranty period shall void all warranty.
KARL STORZ gives no warranty on the correct functioning of units or instruments which have been repaired by unauthorized third parties.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us: all repair work must employ original KARL STORZ parts. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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For rigid replacement telescopes you receive a warranty of 3 years, for fiberscopes and instruments 2 years

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**Important information
for users of KARL STORZ
devices**

**1 Important information
for users of KARL STORZ
devices**

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Thank you very much for your expression of trust in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the UNIDRIVE® S III ARTHRO motor system. The motor system can be used with the following application instruments:

- DrillCut-X® ARTHRO handpiece 28200 DX
- Standard shaver handpiece (8K) 28 721035 or 28 721045 (without head buttons)
- Hybrid shaver handpiece (7K) 28 721050
- High-speed shaver handpiece (15K) 28 721030
- Small joint shaver handpiece (6K) 28 721040
- Multifunction handpiece 28 721036 for use with optional attachments 28 721031/37/38/32/33/34

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Users of the various application instruments should familiarize themselves with the operating instructions included with the respective devices.

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Controls, displays, connectors, and their uses

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3 Controls, displays, connectors, and their uses

- ① Power switch for turning the unit on and off
- ② TFT touch screen
- ③ Connection socket for four-pedal footswitch
- ④ Connection socket for dual-pedal or triple-pedal footswitch
- ⑤ Connection socket for multifunction handpiece
- ⑥ Connection socket for shaver
- ⑦ USB port (for maintenance only)
- ⑧ Potential equalization connection
- ⑨ Power socket
- ⑩ Fuse holder
- ⑪ ® connection sockets*
- ⑫ Dual-pedal footswitch
- ⑬ Cover switch for activating the function which enables the speed to be reduced/increased using the footswitch.

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- ① *** NOTE:** The KARL STORZ SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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1	Important information for users of KARL STORZ devices.....	III
2	Images of the equipment.....	IV
3	Controls, displays, connectors, and their uses.....	V
4	Symbols employed.....	VI
5	General information.....	4
5.1	Description of the system.....	4
5.1.1	Combination of UNIDRIVE® S III ARTHRO and ARTHROPUMP® Power.....	5
5.2	Special features.....	6
6	Safety instructions.....	7
6.1	Explanation of warnings and cautions.....	7
6.2	Normal use.....	13
6.3	Intended use – motor control unit.....	14
6.4	Intended use, indication/ contraindications – accessories for the UNIDRIVE® S III ARTHRO.....	14
6.4.1	Drill/hollow-core drill/burr.....	14
6.4.2	Drill wire/Kirschner wire.....	15
6.4.3	Multifunction handpiece/ drill handpiece/drive incl. attachments.....	16
6.4.4	Shaver handpiece.....	16
6.4.5	Shaver blades.....	17
6.4.6	Saw attachments.....	18
6.5	User qualifications.....	18
6.6	User profile (specialist physician).....	19
6.7	User profile (assistant).....	19
6.8	Patient profile.....	19
6.9	Affected body part or tissue type.....	20
6.10	Use.....	20
6.11	Safety precautions at the site of installation.....	21
6.12	Safety precautions when operating the system.....	21
7	Operating instructions.....	23
7.1	Unpacking.....	23
7.2	Basic equipment.....	23
7.3	Installing and connecting the unit.....	24
7.4	KARL STORZ SCB.....	25
7.5	Connecting the footswitch.....	25
7.6	Connecting the multifunction handpiece/shaver handpiece.....	26
7.7	DrillCut-X® ARTHRO handpiece (28200 DX).....	27
7.7.1	Disassembling the valve assembly from the handpiece.....	27
7.7.2	Assembling the valve assembly on the handpiece.....	27
7.7.3	Assembling the blade on the handpiece.....	28

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7.8	Attaching the blade to the shaver handpiece (6K, 8K, 15K)	29
7.9	Attaching the blade to the 7K shaver handpiece 28 7210 50	29
7.10	Connecting the suction tube (6K, 7K, 8K, 15K, 28200 DX)	30
7.11	Controlling suction (6K, 7K, 8K, 15K, 28200 DX)	30
7.12	Operation of the handpieces (DrillCut-X® ARTHRO, 6K, 8K and 15K)	31
7.13	Handpiece button operation (7K shaver handpiece 28 7210 50)	32
7.14	Small joint shaver handpiece (6K)	33
7.15	Shaver handpiece 28721045	33
7.16	Footswitch operation	34
7.16.1	Dual-pedal footswitch	34
7.16.2	Triple-pedal footswitch	35
7.16.3	Four-pedal footswitch	35
7.17	Multifunction handpiece (drill/saw)	36
7.18	Functional check of the shaver blades before renewed application	37
7.19	Initial operation	38
7.20	UNIDRIVE® S III function	38
7.21	Selecting the language	39
7.22	Selecting the motor connection	40
8	Quick overview of screen functions	41
9	Mode: Settings	43
9.1	System information	43
9.2	Language	44
9.3	Date/Time	44
9.4	Service	44
9.5	Error messages	45
10	Maintenance	46
10.1	Changing the fuse	46
10.2	Reprocessing of the device	47
10.2.1	General warnings	47
10.2.2	Manual wipe-down disinfection	47
10.3	Reprocessing of the DrillCut-X® ARTHRO and shaver blades	48
10.3.1	General warnings	48
10.3.2	Accessories	49
10.3.3	Disassembly	50
10.3.4	Preparation for cleaning and disinfection	51
10.3.5	Manual precleaning	51
10.3.6	Manual cleaning (only applies to shaver blades)	52
10.3.7	Manual disinfection (only applies to shaver blades)	52
10.3.8	Machine cleaning and disinfection	53

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10.3.9	Assembly, inspection and care.....	53
10.3.10	Packaging systems.....	55
10.3.11	Sterilization.....	55
10.3.12	Reprocessing limits.....	55
10.4	Maintenance and safety check.....	56
10.4.1	Maintenance.....	56
10.4.2	Safety check.....	56
10.5	Service and repair.....	58
10.6	Disposal.....	58
10.7	Repair program.....	59
10.8	Important information.....	59
10.9	Limitation of liability.....	60
10.10	Manufacturer's warranty.....	60
11	Technical description.....	61
11.1	General.....	61
11.2	UNIDRIVE® S III ARTHRO error table.....	63
11.3	Troubleshooting.....	65
11.4	Technical Data.....	67
11.5	Intended handpieces.....	68
11.6	Standard conformity.....	69
11.7	Directive compliance.....	69
11.8	Technical documentation.....	70
11.9	Block diagram.....	71
12	Spare parts, recommended accessories	72
12.1	Spare parts.....	72
12.2	Recommended accessories.....	73
12.3	Recommended shaver blades and burrs.....	74
12.4	Recommended shaver blades and burrs for small joint handpiece.....	86
12.5	Accessories for shaver handpieces and SL blades.....	89
12.6	Accessories for pin driver.....	89
12.7	Accessories for sagittal saw.....	89
13	Single-use blades and burrs.....	90
13.1	Intended use.....	90
13.2	Use of the single-use products and explanation of the symbols employed.....	90
14	Electromagnetic Compatibility (EMC) Information.....	92
15	Subsidiaries.....	106

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5.1 Description of the system

The UNIDRIVE® S III ARTHRO can be used with:

- DrillCut-X® ARTHRO handpiece 28200 DX
- Standard shaver handpiece (8K) 28721035
- Standard shaver handpiece without head buttons (8K) 28721045
- Hybrid shaver handpiece (7K) 28721050
- High-speed shaver handpiece (15K) 28721030
- Small joint shaver handpiece (6K) 28721040
- Multifunction handpiece 28721036
- Dual-pedal footswitch 20016831 or Three-pedal footswitch 20012832

It is used in conjunction with a shaver blade or burr (see Accessories page 74-88) to perform surgical interventions under endoscopic control in the field of arthroscopy.

With the corresponding blade, the system is designed to provide controlled cutting, shaving, and abrading of bone and tissue during arthroscopic surgical procedures of the knee, shoulder, ankle, elbow, wrist, hip, and TMJ.

The ergonomically designed and lightweight shaver handpiece has a suction channel for connection to a suction system. Using a rocker lever, the suction channel can be opened or closed.

The easily visible LCD/touch screen on the control unit of the UNIDRIVE® S III ARTHRO provides information on the most important operating parameters. The touch screen provides control of the operating parameters.

With the optionally available multifunction handpiece and its accessories, such as the drill attachment and the saw attachment, the UNIDRIVE® S III ARTHRO becomes a multifunctional instrument that can be used in a wide range of areas in orthopedics.

The control unit automatically recognizes the shaver handpieces or the multifunction handpiece as soon as they are connected.

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5.1.1 Combination of UNIDRIVE® S III ARTHRO and ARTHROPUMP® Power

Either a dual-pedal footswitch or the combined suction and irrigation pump ARTHROPUMP® Power 28340720 can be connected to the right, lower connection socket ④ (see page IV).

If an ARTHROPUMP® Power is connected and a shaver is connected to position ⑥ on the UNIDRIVE® S III ARTHRO and turned on, the shaver sends the signal 'Shaving in progress' to the ARTHROPUMP® Power (increased suction). The tube squeeze valve of the ARTHROPUMP® Power then switches from the function 'Sheath' to the function 'Shaver suction'.

Please observe the instruction manual for the ARTHROPUMP® Power 96156010DER.

The four-pedal footswitch 20016834 is used if the suction and irrigation pump 'ARTHROPUMP® Power' is connected to the socket ④ on the UNIDRIVE® S III ARTHRO via the connection cable 28340770 (see instruction manual for the ARTHROPUMP® Power).

The two left pedals (gray) perform the same function as with the dual-pedal footswitch (counterclockwise/clockwise rotation).

The blue pedal triggers the lavage function. When the inflow is closed, pressing the pedal activates the suction function.

The red pedal triggers the hemostasis mode (temporary increase in pressure). After a predetermined period of time x has elapsed, but after no more than 2 minutes, the pump returns to normal mode. Pressing the red pedal again during hemostasis mode deactivates this mode.

NOTE: Additional equipment connected to medical electrical equipment must be able to prove its compliance with the respective IEC or ISO standards. Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements.

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5.2 Special features

- 2 operating modes of the shaver handpieces: 'Normal' and 'Oscillation'
- Speed range of the small joint shaver handpiece (6K): 1,000 rpm to 6,000 rpm *
- Speed range of the hybrid shaver handpiece (7K): 1,000 rpm to 7,000 rpm *
- Speed range of the standard shaver handpiece (8K) and DrillCut-X® ARTHRO handpiece: 1,000 rpm to 8,000 rpm *
- Speed range of the high-speed shaver handpiece (15K): 1,000 rpm to 15,000 rpm *
- Speed range of the oscillation mode of the shaver handpieces: 500 rpm to 3,000 rpm
- The speed and operating mode of the shaver handpieces can be controlled by either the touch screen, the handpiece (except 28 7210 40 und 28 7210 45) or the footswitch.
- Fast and simple blade change using a precise, quick-release coupling.
- UNIDRIVE® S III ARTHRO is equipped with the KARL STORZ Communication Bus (SCB).
- The control unit automatically recognizes the shaver handpiece and the multifunction handpiece connected in parallel.

i * **NOTE:** The recommended maximum speed is 4,000 rpm in 'Normal' mode and 1,500 rpm in 'Oscillation' mode.

Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

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Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved.

CAUTION: A Caution indicates that particular maintenance procedures or safety precautions must be followed to avoid possible damage to the device.

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WARNING: Read this instruction manual thoroughly before operating the system. Read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

WARNING: Do not operate the handpiece unless the blade/burr is visible.

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Safety instructions

Warnings and cautions

-  **WARNING:** Continuous checks must be performed both before and during use to establish whether the required channel and/or mode is set.
-  **WARNING:** The UNIDRIVE® S III ARTHRO may only be used by persons with an appropriate medical qualification and who are acquainted with the technique.
-  **WARNING:** When replacing fuses, use only fuses of the same rating.
-  **WARNING:** The device may only be operated at the voltage stated on the rating plate.
-  **WARNING:** If the motor does not stop as required, switch the unit off.
-  **WARNING:** Store/keep the instrument at a safe distance from the patient when not in use.
-  **WARNING:** The UNIDRIVE® S III ARTHRO must only be used in medical facilities that conform to applicable national, state and local electrical codes.
-  **WARNING:** Keep out of reach of patients.
-  **WARNING:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if:
 - they are identified as such in the respective instruction manuals or
 - the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or Para. 16 of the 3rd Edition of IEC 60601-1).
-  **WARNING:** Connect the device to a properly grounded 'Hospital grade' or 'Hospital use' electrical outlet. Check grounds regularly. Routinely inspect the power cord and plug and do not use if the plug or the cord is damaged.

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Safety instructions

Warnings and cautions

-  **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.
-  **WARNING:** Do not use in the presence of flammable anesthetics. There is a risk of explosion.
-  **WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.
-  **WARNING:** Risk of injury: Incorrectly assembled and damaged instruments can lead to burns or other injuries to the patient or user. Instruments and all of the accessories used in combination with them must be checked and if necessary replaced immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.
-  **WARNING:** Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

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Safety instructions

Warnings and cautions

WARNING: Risk of injury: Incorrect application of medical instruments or bending of these instruments by applying excessive side load force poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application

WARNING: Risk of injury and risk of damaging the products: When using applied parts outside of the field of vision there is a risk of unintentionally injuring tissue and damaging accessories. Always keep the applied parts of the active electrode and the laser and other instruments which transmit energy pointed at the target in the field of vision during use.

WARNING: Risk of infection: Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

NOTE: The 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions can be downloaded or requested at www.karlstorz.com.

WARNING: When using HF surgical devices, the shaver handle and shaver blade must be put out of reach of the patient.

WARNING: Always unplug the unit before performing any maintenance, including cleaning.

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Safety instructions

Warnings and cautions

WARNING: To avoid the risk of electrical shock, do not open the device. Refer servicing to the manufacturer or qualified personnel that are authorized by the manufacturer. Removal of the cover by non-authorized personnel will void the warranty.

CAUTION: The UNIDRIVE® S III ARTHRO may only be used with accessories, and disposable items which have been designed for use by KARL STORZ.

CAUTION: Do not store liquids on or above the device. Do not allow fluids to get into the unit.

CAUTION: If the motor becomes blocked for any length of time, it can lead to overheating and cause damage to the motor. If the blockage cannot be disengaged using oscillating operation (if possible), the motor must be shut off and the blockage must be removed manually.

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- i NOTE:** The USB port fitted with a cover is intended for use by Service.
- i NOTE:** A defective sterile glove can result in contamination if it comes into contact with sterile goods and can make replacement necessary.
- i NOTE:** Before using the unit, a check must be performed to establish whether the touch screen and the unit display are functioning correctly.
- i NOTE:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.
- i NOTE:** For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the reprocessing of medical devices. Current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.
- i NOTE:** If the unit is not functioning as required, please consult the display or contact the KARL STORZ Service.
- i NOTE:** Limits can be adjusted when working with RUI (= the SCB software). When working with RUI, please check the parameters before and after use.
- i NOTE:** This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your distributor for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Before using the UNIDRIVE® S III ARTHRO motor system on a patient, it is imperative that the instructions in this manual are understood and are followed carefully.

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6.2 Normal use

The UNIDRIVE® S III ARTHRO motor system is universally applicable in conjunction with the respective, suitable application instrument in the ORTHO field of application. The detailed „Intended use“ of the application instrument can be found in the instructions for use provided for the respective application instrument.

The UNIDRIVE® S III ARTHRO motor system may only be used with accessories, wearing parts and disposable products that are designated by KARL STORZ as suitable for the unit.

The UNIDRIVE® S III ARTHRO motor system is used in conjunction with a shaver blade or burr to perform surgical interventions under endoscopic control in the field of arthroscopy.

With the corresponding blade, the system is designed to provide controlled cutting, shaving, and abrading of bone and tissue during arthroscopic interventions on the knee, shoulder, ankle, elbow, wrist, hip, and TMJ.

The ergonomically designed and lightweight handpiece has a suction channel for connection to a suction system. Using a rocker lever, the suction channel can be opened or closed.

With the multifunction handpiece and its accessories, such as the drill attachment and the saw attachment, the UNIDRIVE® S III ARTHRO becomes a multifunctional instrument that can be used in a wide range of areas in orthopedics.

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When combining the UNIDRIVE® S III ARTHRO and appropriate pump as a system, the instructions for use provided for the respective combination device must also be observed.

Use of the UNIDRIVE® S III ARTHRO in fields other than those indicated above is not permitted.

The UNIDRIVE® S III ARTHRO may only be used with accessories, replacement parts, and disposable items which have been designated as suitable for the unit by KARL STORZ.

Unauthorized conversions or modifications to the system are not permitted.

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6.3 Intended use – motor control unit

When used in arthroscopy/orthopedics, the motor control unit serves to regulate and control motors or handpieces and accessories during orthopedic, arthroscopic and minimally invasive interventions for the excision, removal, ablation and transection of soft tissue and/or bones.

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6.4 Intended use, indication/ contraindications – accessories for the UNIDRIVE® S III ARTHRO

6.4.1 Drill/hollow-core drill/burr

6.4.1.1 Intended use

Drills and burs are used for drilling, milling and separating bone during brief open and arthroscopic interventions in the musculoskeletal system.

6.4.1.2 Indications

- Epi- and metaphyseal bone tunnel
- Bone harvesting for cancellous bone grafts and cartilage-bone transplants

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6.4.1.3 Contraindications

- Contact with the central nervous system
- The accessories for the motor system for musculoskeletal applications must not be used if, in the opinion of the responsible physician, the patient's health is endangered through their use, e.g., due to the patient's general condition.
- Further contraindications relating directly to the product group are not known.

6.4.2 Drill wire/Kirschner wire**6.4.2.1 Intended use**

Wires are used for the fixation of bone fragments and the accurate introduction of cannulated instruments and implants during brief open and arthroscopic interventions in the musculoskeletal system.

6.4.2.2 Indications

- Guidance of cannulated instruments for the creation of epi- and metaphyseal bone canals
- Guidance of cannulated implants
- Reduction and fixation of metaphyseal fractures
- Reduction and fixation of diaphyseal fractures and luxation of hand and foot bones
- Temporary arthrodesis of small joints
- Temporary intraoperative fixation of fracture fragments

6.4.2.3 Contraindications

- Contact with the central nervous system
- The accessories for the motor system for musculoskeletal applications must not be used if, in the opinion of the responsible physician, the patient's health is endangered through their use, e.g., due to the patient's general condition.
- Further contraindications relating directly to the product group are not known.

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**6.4.3 Multifunction handpiece/
drill handpiece/drive incl.
attachments**

6.4.3.1 Intended use

The multifunction handpiece serves to hold and drive various attachments used for drilling and cutting tissue during therapeutic interventions of the musculoskeletal system.

6.4.3.2 Indications

The multifunction handpiece/drill handpiece and corresponding accessories are suitable for all therapeutic interventions which require the energy/media supplied by the device for successful operation.

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6.4.4 Shaver handpiece

6.4.4.1 Intended use

The shaver handpieces serve to hold and drive shaver blades designated by KARL STORZ for the ablation and removal of soft tissue and/or bone during brief arthroscopic interventions.

6.4.4.2 Indications

The shaver handpiece and corresponding accessories are suitable for all diagnostic and therapeutic interventions which require the energy/media supplied by the device for successful operation.

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6.4.4.3 Contraindications

- Contact with the central nervous system.
- The accessories for the motor system for musculoskeletal applications must not be used if, in the opinion of the responsible physician, the patient's health is endangered through their use, e.g., due to the patient's general condition.
- Contraindications relating directly to the product group are not known.

6.4.5 Shaver blades**6.4.5.1 Intended use**

Shaver blades serve to ablate and remove soft and hard tissue in synovial joints during brief arthroscopic interventions.

6.4.5.2 Indications

- Arthroscopic debridement in the knee joint
- Arthroscopic debridement in the shoulder joint
- Arthroscopic debridement in the hip joint
- Arthroscopic debridement in the ankle joint
- Arthroscopic debridement in the wrist joint
- Arthroscopic debridement in the elbow joint
- Arthroscopic bone resection with impingement

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6.4.6 Saw attachments

6.4.6.1 Intended use

Together with the appropriate saw attachment and drive, saw attachments serve to separate bone in the musculoskeletal system during brief surgical interventions.

6.4.6.2 Indications

- Tibial and patellar bone block harvesting for bone-tendon-bone transplants
- Osteotomies, epi-, meta- and diaphyseal bone

6.5 User qualifications

The UNIDRIVE® S III ARTHRO motor system in conjunction with an application instrument may only be used by physicians and/or medical assistants who have been instructed in the use of this device/instrument.

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6.6 User profile (specialist physician)

- Recognized medical skills in the area of arthroscopy (specialist physician)
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

6.7 User profile (assistant)

- Completed training as medical assistant
- Adequate language skills in one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

6.8 Patient profile

The use of the UNIDRIVE® S III ARTHRO is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health with regard to the endoscopic intervention and declare it adequate.

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6.9 Affected body part or tissue type

The area of application is limited to certain joints (feet, knees, hips, shoulders, elbows, hands). The joints are accessed via surgically invasive openings. The area of application must be in a condition suitable for the endoscopic intervention.

6.10 Use

The device may only be used in operating rooms in hospitals and medical practices in the ambient conditions given in the technical data (page 67).

- Not on intensive care stations.
- Not in mobile facilities (ambulances, aircraft, ships or trains).
- Position it on a vibration-free, fixed base.
- Visual distance from the physician to the device when operated: max. 2 m.
- Visual distance for operating the device: 30-70 cm.
- Viewing angle to the front of the device: max. $\pm 45^\circ$.
- The use of the device requires the division of work.
- Frequency of use up to several times a day.
- Length of use from a few minutes to several hours a day.
- Can be moved only by installing it on an equipment trolley.
- Can be swiveled only by installing it on ceiling supply units.
- The patient must not necessarily be under general anesthesia during the intervention (local anesthesia).
- May be used with other devices (light source, shaver, HF surgery, etc.).
- To be operated in ambient light between 100 lux and 1500 lux.
- Ambient temperature in immediate proximity of the device during operation: max. 40°C.
- Sound pressure from background noise < 60 dBA.
- The device is placed outside the sterile area.

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
**6.11 Safety precautions at the site of installation**

The UNIDRIVE® S III ARTHRO motor system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The UNIDRIVE® S III ARTHRO is equipped with a connector for attaching a ground wire. Connect this in accordance with the state, local and national electrical codes.

The UNIDRIVE® S III ARTHRO motor system is not intended for use in hazardous areas. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

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6.12 Safety precautions when operating the system

WARNING: The system must be tested for operational safety and proper condition prior to each surgical procedure.

During surgery with the UNIDRIVE® S III ARTHRO, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

CE marked equipment has been tested and is found to comply with the EMC limits of IEC 60601-1-2 CISPR 11 Class B.

Please observe the Electromagnetic Compatibility information in Chapter 14 (p. 92-104).

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WARNING: When using HF surgery, the shaver handpiece and shaver blade must not come in contact with the patient.

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7 Operating instructions

7.1 Unpacking

Carefully remove the unit and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

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If possible, retain the original packing materials for later use; these can be used if the unit has to be transported.

7.2 Basic equipment

- 1 UNIDRIVE® S III ARTHRO 28 7230 20-1
- 1 Power cord 400 A
- 1 SCB connecting cable, length 100 cm, 20 0901 70
- 1 Instruction manual 96156018 DER

ⓘ NOTE: Optional accessories are described in section 'Recommended accessories'.

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
NOTE: The UNIDRIVE® S III ARTHRO including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

! WARNING: The UNIDRIVE® S III ARTHRO is not intended for use in hazardous areas. Do not operate the system within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
Place the UNIDRIVE® S III ARTHRO on a level surface.

The UNIDRIVE® S III ARTHRO is equipped with a connector ⑧ for attaching a ground wire. The unit's grounding should be installed by a qualified electrician.

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Before connecting the UNIDRIVE® S III ARTHRO to the power supply line, verify that the line voltage stated on the unit's identification plate conforms with that of the power supply lines to be used.

Connect the power, insert the power cord end into the receptacle ⑩ of the unit.

! WARNUNG: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

! WARNUNG: Only use the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

! WARNUNG: The device must be set up so as to ensure that the power cord can simply be pulled at any time.

! WARNUNG: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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
7.4 KARL STORZ SCB

① NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector is equipped with a protection mechanism. Pull the protection mechanism of the SCB connector and insert the plug into one of the SCB sockets ⑪. Connect the other end of the cable to a KARL STORZ SCB control unit (KARL STORZ Communication Bus) or other SCB devices (see instruction manual KARL STORZ SCB control NEO System).

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7.5 Connecting the footswitch

To use the footswitch ⑫, insert the connecting cord of the footswitch into the socket for the footswitch (③ or ④) on the front of the UNIDRIVE® S III ARTHRO unit. Align the red dot on the connecting cord with the red dot on the footswitch socket.

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
7.6 Connecting the multifunction handpiece/shaver handpiece

Prepare the multifunction handpiece/shaver handpiece as described in the instruction manual.

NOTE: Prior to turning the device on, the user must connect a multifunction handpiece and/or shaver handpiece to the device. This is checked by the software and, if nothing is connected, a corresponding error message is shown in the top image margin.

NOTE: Please follow the instruction manual enclosed with the multifunction handpiece/shaver handpiece!

Connect the connecting cord from the multifunction handpiece/shaver handpiece to the connection socket ⑤ (multifunction handpiece) and/or ⑥ (shaver).

Align the red point on the connecting cord with the red point on the connection socket to guarantee a proper connection.

NOTE: If no 'motor' is connected when the unit is switched on, the message: 'E12: NO MOTOR' flashes in the status line of the display.

NOTE: To unplug the connecting cable, pull out the plug — do not pull directly on the cable!

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7.7 DrillCut-X® ARTHRO handpiece (28200 DX)

7.7.1 Disassembling the valve assembly from the handpiece

① **NOTE:** The valve assembly comprises the valve lock ① and valve body ② (see Fig. A).

1. Turn the valve lock ① counterclockwise to 11 o'clock (see Fig. B).
2. Remove the valve assembly axially from the handpiece by pulling toward distal (see Fig. C).
3. Dismantling the valve assembly:
Turn the valve lock ① counterclockwise to 9 o'clock and remove from the valve body ② by pulling (see Fig. D and E).

7.7.2 Assembling the valve assembly on the handpiece

1. Place the valve lock ① on the valve body ② (90° to the left toward the valve lever) (see Fig. E and F).
2. Continue to turn the valve lock ① clockwise to 11 o'clock until the valve lock noticeably clicks into place (see Fig. G) with the locking ball (see arrow in Fig. F).
3. Push the valve assembly toward proximal on the axis of the handpiece as far as it will go, press on firmly (see Fig. C).
4. Lock the valve assembly on the handpiece:
Turn the valve lock ① clockwise to 12 o'clock (the valve lock and valve body are aligned, see Fig. H).
5. Pull axially toward distal to check that the valve assembly is securely in place.

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



7.7.3 Assembling the blade on the handpiece

1. Pull the locking sleeve on the handpiece back toward proximal and hold (in this position) (see Fig. I).
2. Push the blade in the desired position (see note) into the handpiece (blade mount) until it locks into place (stop) (see Fig. I und J).

NOTE: The 'blade window' (see Fig. K) can be fixed radially in 4 positions (90°, 180°, 270° and 360° to the valve lever).

3. Let go of the locking sleeve (blade mount) and then pull axially on the blade to ensure it is securely in place (see Fig. J/L).

CAUTION: Danger of damage to the valve assembly!
To ensure the proper functioning of the valve assembly (valve lock, valve body) in the long term, it must be handled with care.

WARNING: When the valve assembly is not assembled, ensure that the handpiece cannot be activated.

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NOTE: Removing the blade from the shaver handpiece: In a deactivated state, pull back the locking sleeve on the handpiece and pull the blade out of the handpiece (see Fig. J).

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7.8 Attaching the blade to the shaver handpiece (6K, 8K, 15K)

NOTE: Insert the blade either with direction setter, and thus with the blade opening facing upwards, or rotated through 180°, with the opening facing downwards.

1. Insert the blade with the direction setter into one of the two notches on the lock.
2. Turn the lock on the lever counterclockwise and hold it.
3. Slide the blade **as far as it will go** (see arrow) into the lock.
4. Finally, let go of the lever and make sure that the blade is securely attached by pulling it.

NOTE: During use, ensure that the locking mechanism does not become unlocked.

NOTE: If disposable products are used, follow the instructions on pages 90ff.

NOTE: Removing the blade from the shaver handpiece: In a deactivated state, turn the lock on the lever counterclockwise and hold in this position. Pull the blade out of the handpiece. During use, ensure that the locking mechanism does not become unlocked.

NOTE: The shavers stop automatically after a continuous operating period of approx. 10 minutes. Releasing and pressing the footswitch pedal or button on the handpiece again restarts the shaver.

7.9 Attaching the blade to the 7K shaver handpiece 28721050

Press the button on the side of the handpiece collar and insert the selected shaver blade into the receptacle on the handpiece until the shaver blade clicks into place.

NOTE: Be sure to align the notch on the blade with the corresponding notch on the handpiece.

NOTE: During use, ensure that the locking mechanism does not become unlocked.

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
ⓘ NOTE: Removing the blade from the 7K shaver handpiece 28721050: In a deactivated state, release the shaver blade by pressing the button on the side of the shaft of the handpiece and removing the shaver blade.

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
7.10 Connecting the suction tube (6K, 7K, 8K, 15K, 28200 DX)

Attach the suction tube to the connection on the handpiece and connect the tubing to a suction pump or wall suction. Suction is controlled using the rocker lever. If the rocker lever is at the rear limit stop, minimum suction is provided, at the front limit stop, maximum suction is provided.

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
7.11 Controlling suction (6K, 7K, 8K, 15K, 28200 DX)

Suction is controlled using the rocker lever (see arrow). If the rocker lever is at the rear limit stop, the suction channel is closed; at the front limit stop, the suction channel is open (maximum suction).

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⚠ CAUTION: During shaver activation the handpiece needs a certain amount of flow for the cooling of the handpiece. Make sure that the irrigation/suction system is properly set up and that tubes are not occluded or kinked.

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
⚠ CAUTION: The minimum suction position restricts or eliminates fluid movement. Operating the handpiece with minimal fluid movement can result in excess heat retention and damage to the device.

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
7.12 Operation of the handpieces (DrillCut-X® ARTHRO, 6K, 8K and 15K)

The UNIDRIVE® S III ARTHRO system features handpieces which can be operated via the LCD/touch screen of the UNIDRIVE® S III ARTHRO, the buttons of the handpiece (except 28 721040 and 28 721045) and the footswitch. The handpieces 28200 DX, 28 721030 and 28 721035 have three buttons with functions which are identical on all three handpiece types. The button functions are explained in the following:

Button 1/III (right): Oscillation function or »+«

- Non-running handpiece: Pressing this button switches to Oscillation mode (change in direction). Each time the button is pressed again, the blade oscillates for 1 second.
- Running handpiece: Pressing the button 1/III increases (+) the speed.

Button 2/II (left): Change direction of rotation function (forward/reverse) or »—«

- Non-running handpiece: Pressing the button 2/II changes the mode of rotation direction from forward to reverse.
- Running handpiece: Pressing the button 2/II reduces the speed.

Button 3/I (center): Start/Stop function

This button starts and stops the handpiece.

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7.13 Handpiece button operation (7K shaver handpiece 28 7210 50)

The UNIDRIVE® S III ARTHRO system features a fourth shaver handpiece (7K), which can also be controlled via the touch screen/LCD, handpiece buttons, or footswitch. The 7K handpiece is equipped with two sets of buttons. The first set consists of three buttons as depicted on the left. These buttons are located at the top of the handpiece. The second set consists of a single button located at the bottom of the handpiece. This button is used to start and stop the handpiece.

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Change direction of rotation function (forward/reverse) or »«

- With the handpiece stopped, pressing this button will change the direction of the blade rotation from forward to reverse.
- When the handpiece is running, pressing this button will decrease the speed of the blade.

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
Oscillation function or »+«

- When the handpiece is stopped, pressing this button will switch to oscillation mode (change of direction). Each time the button is pressed again, the blade oscillates for 1 second.
- When the handpiece is running, pressing this button will increase the speed of the blade.

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
7.14 Small joint shaver handpiece (6K)

The small joint shaver handpiece is operated using either the LCD/touch screen of the UNIDRIVE® S III ARTHRO system or the footswitch. It only has one button, with the following function:

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7.15 Shaver handpiece 28721045

The shaver handpiece 28721045 is operated using either the LCD/touch screen of the UNIDRIVE® S III ARTHRO system or the footswitch.

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7.16 Footswitch operation

The shaver handpieces can be controlled using a footswitch.

NOTE: It is recommended to keep a footswitch ready as a backup solution in the event of failure of the control buttons.

NOTE: The shaver and the multifunction handpiece stop automatically after a continuous operating period of approx. 10 minutes (shaver)/after approx. 5 minutes (multifunction handpiece). Releasing and pressing the footswitch pedal or button / trigger on the handpiece again restarts the shaver or multifunction handpiece.

7.16.1 Dual-pedal footswitch

The dual-pedal footswitch 20016831 has three 'controls': a left pedal, a right pedal and a cover switch.

Below is a definition of the footswitch functionality:

Left pedal:

Reverse (Left) < Counterclockwise

Right pedal:

Forward (Right) > Clockwise

Left and right pedal simultaneously:

Oscillate 

Cover switch (button on top of the footswitch):

Pressing the cover switch once switches the footswitch into speed adjust mode; the foot pedals can now adjust the speed up and down.

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By pressing the cover switch again, the device returns to normal operating mode with the new speed setting.

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
7.16.2 Triple-pedal footswitch

The triple-pedal footswitch 20012832 has three distinct pedals for operating mode selection and a switch for speed control.

Below is a definition of the triple-pedal footswitch functionality:

Left pedal:

Reverse (Left) < Counterclockwise

Middle pedal:

Oscillate (○)

Right pedal:

Forward (Right) > Clockwise

Button on top of footswitch:

Pressing the button on the top of the footswitch once switches the footswitch into speed adjust mode; the foot pedals can now adjust the speed up and down.

Pressing the right foot pedal increases the speed (by the programmed rate) and pressing the left foot pedal decreases the speed (by the programmed rate).

By pressing the top button again, the device returns to normal operating mode with the new speed setting.

7.16.3 Four-pedal footswitch

The four-pedal footswitch 20016834 is used if the suction and irrigation pump 'ARTHROPUMP® Power' 28340720-1 is connected to the socket ④ on the UNIDRIVE® S III ARTHRO via the connection cable 28340770 (see instruction manual for the ARTHROPUMP® Power).

The two left pedals (gray) perform the same function as the dual-pedal footswitch (counterclockwise/clockwise rotation). The blue pedal triggers the lavage function. When the inflow is closed, pressing the pedal activates the suction function.

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
28721033

28721034

7.17 Multifunction handpiece (drill/saw)

The UNIDRIVE® S III ARTHRO control unit is capable of driving a multifunction handpiece. This handpiece may be attached to drill chucks, sagittal saws, etc. A list of available accessories is shown in a separate section of this manual.

Each accessory is connected to the handpiece by attaching the accessory and handpiece together until solid contact is made. Press the rear of the button on the top 17 of the handpiece to remove the accessory.

The motor is cannulated and accepts pins up to 0.125 inches (3.2 mm) in diameter.

The UNIDRIVE® S III ARTHRO allows the speed of the handpiece to be controlled by moving the trigger 16. The speed (rpm) of the handpiece is proportional to the position of the trigger. Pushing the trigger all the way in toward the handpiece provides the maximum speed.

NOTE: The multifunction handpiece stops automatically after a continuous operating period of approx. 5 minutes. Releasing and pressing the footswitch pedal or trigger on the handpiece again restarts the multifunction handpiece.

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



7.18 Functional check of the shaver blades before renewed application

WARNING: Risk of injury to the patient and/or user: Before every renewed application of the shaver blades, it is compulsory to carry out a functional check. This includes the following tests:

- The inner blade should be free to rotate easily within the outer blade and must not show any run-in marks (negative example, see photos 1 and 2).
- The inner blade must lock securely into place in its front end position, in order to ensure that it is not unintentionally lost.
- The distal cutting edges must be free from defects (negative example, see photo 3).
- The outer blade of the burr must not show any points of abrasion (negative example, see photos 3 and 4).

Always operate the shaver handpiece/shaver blades with running suction

WARNING: Risk of injury to patient and damage to property:

It is strongly recommended never to operate the shaver handpieces/shaver blades without running suction. Suction of the irrigation fluid serves to cool the running surfaces and the shaver motor and prevents injuries/damage caused by overheating and mechanical 'seizure'. To this end, the rocker lever used for setting the suction power must at least be in a slightly open position (see photo 5).

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
7.19 Initial operation

Switch the UNIDRIVE® S III ARTHRO on at the power supply ① (power switch lights up).

! NOTE: Prior to turning the device on, the user must connect a multifunction handpiece and/or shaver handpiece to the unit. This is checked by the software and if necessary (if nothing is connected) indicated by a corresponding error message in the top image margin.

The screen shown on the left appears for the time it takes the application to start.

! NOTE: In addition, the following pages describe the simplest application case as an example, as well as the adjustment of the speed. A single signal tone can be heard during the switch-on procedure.

7.20 UNIDRIVE® S III function

A start screen is displayed once the unit has started. The function only becomes active once the multifunction handpiece and/or shaver handpiece has/have been connected to the unit.

! CAUTION: In order to avoid damaging the delicate touch screen interface, do not use any objects to touch the touch screen. The screen displayed shows the start screen with a shaver handpiece connected.

! NOTE: Error messages are shown in text in the status line of the operation screen if necessary (flashing text).

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7.21 Selecting the language


It may be required to define a different language for the operation screen since the device is set to English as the default in the factory.

This is done by pressing the 'menu' symbol (wrench) in the bottom left.

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Activate the 'Language' button (2nd button from the top).

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Select the desired language.

Select one of the ten available languages by pressing the respective language button. The languages currently available are displayed over two screens.

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7.22 Selecting the motor connection

The handpieces described on page 4 (which must be operated within various speed ranges) can be connected up to the control unit. To do so, the required instrument must be connected to the corresponding motor connection socket (connect multifunction handpiece up to upper left socket ⑤, connect shaver handpiece up to upper right socket ⑥).

Having detected the handpiece automatically, the UNIDRIVE® S III ARTHRO loads the standard parameters and displays them on the screen. The symbol for the respective instrument is displayed in the bottom right of the screen.

It is only possible to switch between the individual motor connection sockets when both channels are occupied.

The switching is performed via the forward/backward lever on the multifunction handpiece. In the 'OFF' setting, the shaver handpiece is active, otherwise the multifunction handpiece is active.

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
**8 Quick overview of screen
functions**

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❶ Name of the currently active instrument (multifunction handpiece/shaver handpiece).

Error message E12 appears ('No motor') if neither a multifunction handpiece nor a shaver handpiece is connected.

If both a shaver handpiece and a multifunction handpiece are connected to the unit, the multifunction handpiece is automatically active unless the switch on the lower end of the multifunction handpiece is set to OFF. Then the shaver handpiece is active.

❷ Currently set speed limit, i.e., the maximum number of revolutions currently possible.

The speed limit of the currently active instrument is automatically recognized and displayed by the unit.

❸ Increase/reduction of speed limit in increments dependent on instrument, e.g., 100 rpm for multifunction handpiece.

❹ Display of speed range possible for the currently active instrument.

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Quick overview of screen functions

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5 Footswitch symbols

If the multifunction handpiece is active, the + or - sign appears on the footswitch symbols.

This means that the speed limit 2 on the connected footswitch can be increased/reduced. (For shaver handpieces the cover switch on the dual-pedal footswitch must be pressed first.)

6 Symbols for the forward rotation/oscillation/reverse rotation modes (see also pages 31/32)

Currently active mode = orange.

The 'Oscillation' symbol only appears if a shaver handpiece is currently activated (e.g., when the name of a shaver handpiece is in the header 1).

It is possible to switch between the forward and reverse rotation modes

- via the forward/reverse head buttons of the shaver handpiece
- via the switch at the lower end of the multifunction handpiece (FWD – forward, REV – reverse)

The shaver handpiece can be switched to the oscillation mode via the 'Oscillation' head button.

7 This button can be used to enter the 'Settings' menu. This area allows you to change the screen language (e.g., from English to German) among other functions.

8 The currently active instrument (multifunction handpiece/shaver handpiece) is displayed here in orange. If both instruments are connected to the unit, the multifunction handpiece is automatically active unless the switch on the lower end of the multifunction handpiece is set to OFF.

9 Activate/Deactivate the head buttons of the shaver handpiece (orange = active). This button can thus also be used to turn off a running shaver handpiece.

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9 Mode: Settings

9.1 System information

To call up the menu functions, press the 'Menu' button (wrench) in the bottom left.

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This opens the 'MENU' screen (image on the left). Press the 'SYSTEM INFORMATION' button to obtain detailed information on your UNIDRIVE® S III ARTHRO.

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The device number (serial number) is displayed at the top.


The operating hour counter and the software version statuses are displayed in the lower section.

Touch the tick icon (✓) at the bottom right several times to access the operating menu.

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Siehe Seite 39.

9.2 Language


Here you can set the language to be displayed on the screen.

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9.3 Date/Time

Here you can set the date and time, which appear at the top of the screen.

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
9.4 Service

This menu is reserved for the authorized Service and, for this reason, access to this menu is password-protected.

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The upper section displays the current date and time; the lower section displays a list of all the error messages that have occurred in order of importance.

E01-E10: Hardware error messages which are potentially caused by a unit fault.

E11-E30: User error messages, which are due to incorrect configuration of the unit (e.g. if neither a multifunction handpiece nor shaver handpiece is connected). See Error table on p. 63.

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E 12: NO MOTOR
(No motor/handpiece connected)

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10 Maintenance


10.1 Changing the fuse

Switch the control unit off and disconnect the power.

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
Remove the power fuse holder ⑩ with a screwdriver or another suitable tool.

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! **CAUTION:** Use only fuses with the correct nominal ratings.

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
Reconnect the power.
Perform function testing.

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10.2 Reprocessing of the device

10.2.1 General warnings

! WARNING: Risk of fatal injury from electric shock: Unplug from power supply before every cleaning!

! WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

! WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

! WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

! WARNING: National laws and regulations must be observed.

i NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

10.2.2 Manual wipe-down disinfection

! CAUTION: Damage to the control unit by liquid: Avoid allowing liquid to enter the control unit.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents

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must not be used. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

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• This reprocessing step is part of the validated procedure.

*** NOTE:** Reprocessing instructions for other UNIDRIVE® S III ARTHRO handpieces can be found in the KARL STORZ instruction manual 96156013 DER.

10.3 Reprocessing of the DrillCut-X® ARTHRO and shaver blades

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Maintenance

Reprocessing

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

10.3.2 Accessories

Accessories required for carrying out reprocessing:

- Brushes: 27652, 27650 E
- Irrigation adaptor for shaver blades: 28205 SAA, 28205 SAI

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10.3.3 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible.

1. Pull the locking sleeve on the handpiece back toward proximal and hold (in this position) (see Fig. A).
2. Pull the blade out of the handpiece (blade mount). Release the locking sleeve again (see Fig. B).
3. Remove the suction tube from the suction port on the DrillCut-X® ARTHRO handpiece (see Fig. C).
4. Turn the valve lock ① counterclockwise to 11 o'clock (see Fig. D).
5. Remove the valve assembly axially from the handpiece by pulling toward distal (see Fig. E).
6. Dismantling the valve assembly:
Turn the valve lock ① counterclockwise to 9 o'clock and remove from the valve body ② by pulling (see Fig. F/G).
7. Remove the working insert/inner blade (the inner sheath) of the reusable shaver blade from the outer blade (outer sheath) (see Fig. H).

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

10.3.4 Preparation for cleaning and disinfection

CAUTION: Risk of burns: Dried residue can cause malfunctions and generate heat. The handpiece and the reusable shaver attachments should be subjected to a reprocessing process immediately after use in order to avoid the drying out of tissue remnants and fluids.

WARNING: Reprocessing is only permitted for reprocessable shaver attachments. If you reprocess disposable products, the product safety cannot be guaranteed as a result of material changes. This can result in injury to the patient. Dispose of disposable products after use. Do not reprocess them (also see section 13).

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the handpiece and shaver blades used by wiping down and rinsing, for example. Rinse the shaver blade through, preferably with distilled water (with the cutting opening open) directly after use, in order to prevent the shaver blade from becoming blocked. To do so, activate the suction function. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

10.3.5 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed from the surfaces of the individual components by cleaning them under cold running water with the aid of a brush or sponge (see Fig. a and b).

Brushing the lumina

The lumina and hollow spaces must be precleaned under cold running water with the aid of a suitable brush (see Fig. c and d).

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Removing possible blockages from the valve or blade drive shaft:

Flush the valve body with a syringe.

Rinsing out with a pressure water gun

The lumina and hollow spaces of the medical device must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660; see Fig. e) is ideal for this purpose.

10.3.6 Manual cleaning
(only applies to shaver blades)

The shaver blade must be completely immersed in a cleaning solution. Through disassembling the instrument and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present (see Fig. f). After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

10.3.7 Manual disinfection
(only applies to shaver blades)

The shaver blade must be completely immersed in a disinfectant solution. Through disassembling the instrument and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. Following the exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings and lumina are dried completely with (preferably sterile) compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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10.3.8 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

Connecting

In order to ensure effective machine cleaning and disinfection, the instruments must be connected up to the washer and disinfecter. Connect the DrillCut-X® ARTHRO handpiece up to the washer and disinfecter via a tube and the shaver blades via the irrigation adaptor (see Fig. g).

10.3.9 Assembly, inspection and care

The cleaned and disinfected medical devices must be visually inspected to ensure that they are clean, complete, dry and have not been damaged:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.

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- **Only applies to shaver blades:**
Glide surfaces must be lubricated locally with instrument oil*.
- 1. To do this, apply a small amount of instrument oil to the running surfaces.
- 2. Insert the inner blade into the outer blade.
- 3. Turn the inner blade several times.
- 4. Disassemble the shaver blade again.
- Dismantled medical devices must be assembled. (see section 7.7.2 or 7.7.3)
- Afterwards, a functional check must be carried out.

- ! NOTE:** The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).
- ! NOTE:** The new generation of blades, which can be recognized by the black plastic driver, renders the use of O-rings and washers unnecessary.
- ! NOTE:** Care of the inner blade with instrument oil prior to sterilization increases the service life of the shaver blade.
- ! NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

! WARNING: The DrillCut-X® ARTHRO handpiece must **not** be treated with instrument oil/instrument spray.

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10.3.11 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and following consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization in the fractionated prevacuum procedure

For sterilization of the medical device in its fully assembled state, a fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134°C – 137°C with a minimum exposure time of 4 minutes and a maximum of 18 minutes must be used. This procedure is only suitable for thermostable instruments.

WARNING: *Risk of burns! Parts are hot after steam sterilization. Allow them to cool down!*

WARNING: *KARL STORZ recommends that, in order to assure sterility, all blades/ burs should be sterilized in a disassembled state and then reassembled in a sterile field.*

10.3.12 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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10.4.2 Safety check

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

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Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

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Test for proper functioning

The test for proper functioning must include the following tests at least:

- Test for proper functioning of the power switch
- Footswitch detection
- Handpiece detection
- Test for proper functioning of the handpiece and pump
- Change of direction of rotation
- Reading out system information

Documentation

The safety check and the results must be recorded in the medical device logbook.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during this safety check/repeat inspection, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety checks, please refer to the latest issue of the relevant Service Manual.

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10.5 Service and repair

Repairing defective devices is only to be performed by personnel that have been authorized by KARL STORZ and using original KARL STORZ parts.

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10.6 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

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Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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10.7 Repair program

For devices, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan, which you will then return to KARL STORZ as soon as you receive the repaired unit.

In Germany, you can refer the repairs directly to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

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10.9 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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10.10 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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11 Technical description

11.1 General

In the event of the unit malfunctioning, the troubleshooting list provided in this section will assist the user in finding the possible cause. The troubleshooting list indicates how to find and rectify the various malfunctions.

WARNING: High voltage is present in the unit's interior. The unit must only be opened by an authorized KARL STORZ customer service technician.

WARNING: Maintenance work that is not mentioned in this section must only be performed by technical personnel that are authorized by KARL STORZ.

General information

Error messages are output directly to the screen on the UNIDRIVE® S III ARTHRO. Send the device in for servicing when a hardware error message is displayed (E01-E10). Take appropriate action for user error messages (E12-E30). E.g. E12 'NO MOTOR' – check whether motor/handpiece are properly connected. Detailed information can be obtained from the error table on page 63 (see section 'Service mode' as well).

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**11.2 UNIDRIVE® S III ARTHRO
error table**

Error message	Cause	Remedy
1 = E01: MEMORY ERROR	Parameters could not be saved successfully.	Turn device on/off. Load default values. Or send unit in for repair.
2 = E02: MOTORBOARD DEFECTIVE	Motor control overloaded.	Turn device on/off. Or send unit in for repair.
3 = E03: HARDWARE CONFIGURATION ERROR	Hardware configuration is incorrect.	Send unit in for repair.
4 = E04: SOFTWARE CONFIGURATION ERROR	Software configuration is incorrect.	Send unit in for repair.
5 = E05: –		
6 = E06: –		
7 = E07: –		
8 = E08: –		
9 = E09: MOTOR DEFECTIVE M1	Motor on channel 1 is defective.	Release footswitch and restart. If the error persists, a different handpiece must be used.
10 = E10: MOTOR DEFECTIVE M2	Motor on channel 2 is defective.	See E09
12 = E12: NO MOTOR	No motor/handpiece connected.	Connect motor/handpiece.
13 = E13: TOO HIGH TEMPERATURE MOTOR M1	Motor/handpiece on channel 1 overheated.	Allow the connected handpiece to cool down or operate an additional handpiece on the 2nd channel.
14 = E14: –		
15 = E15: TOO HIGH TEMPERATURE MOTOR M2	Motor/handpiece on channel 2 overheated.	See E13
17 = E17: MOTOR NOT SUPPORTED M1	The motor/handpiece connected on channel 1 is not supported.	Only use handpieces approved by KARL STORZ. A service update may be necessary.
18 = E18: MOTOR NOT SUPPORTED M2	The motor/handpiece connected on channel 2 is not supported.	See E17
19 = E19: MOTOR STOPPED M1	Motor/handpiece on channel 1 blocked in operation.	Release footswitch and restart. If the error persists, clean the handpiece or exchange the attachment.
20 = E20: MOTOR STOPPED M2	Motor/handpiece on channel 2 blocked in operation.	See E19
21 = E21: NO SCB MODULE	Unit is not connected to an active SCB bus.	Insert SCB cable and activate SCB bus.

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Problem

Increased operating noise at the handpiece.

Possible cause

Wear.

Recommended remedy

The intervention can be completed. Send in the handpiece for maintenance.

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11.4 Technical Data

UNIDRIVE® S III ARTHRO	28723020-1
Line voltage	100...240 VAC
Power frequency	50/60 Hz
Line fuse	2 x T 4.0 AH 250 V, IEC 127 format
Power consumption	150...200 W
Dimensions	305 x 165 x 280 mm (W x H x D)
Weight	5.2 kg
Operating mode: periodic duty (in connection with application instrument*)	
ⓘ NOTE: The shaver and the multifunction handpiece stop automatically after a continuous operating period of approx. 10 minutes (shaver)/after approx. 5 minutes (multifunction handpiece). Releasing and pressing the footswitch pedal or button/trigger on the handpiece again restarts the shaver or multifunction handpiece.	
* Shaver handpieces	
Duty cycle:	50 %, t = 10 min (total time 1.5 h)
* Multifunction handpiece	
Duty cycle:	66 %, t = 3 min (total time 10 min)
Operating conditions:	
Temperature	10 °C...30 °C
Humidity	30 % ... 70 %
(non-condensing)	
Atmosph. pressure	700 hPa...1080 hPa
Storage/transport conditions:	
Temperature	-10 °C...60 °C
Relative humidity	20 %...95 %
(non-condensing)	
Atmosph. pressure	500 hPa...1080 hPa
Average lifecycle	
UNIDRIVE® S III ARTHRO	5 years
DrillCut-X® ARTHRO	2 years/max. 300 reprocessing cycles

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6K/7K/8K/15K (oszillierend): 1500 min⁻¹
6K/7K/8K/15K (rotierend): 4000 min⁻¹

**Maximum speed when turned on
Shaver handpiece**

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(for 28723020-1)

According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electrical shock:
Protection class I
- Degree of protection against electric shocks:
Applied part type BF
- Type of moisture protection:
drip-water protection as per IPX 1

Conforms with IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the section 14 (p. 92-104).

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NOTE: The code number after the CE mark indicates the responsible notified body.

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11.8 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

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① NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
11.9 Blockschaltbild

11.9 Block diagram

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**Spare parts,
recommended accessories**

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**12 Spare parts,
recommended accessories**

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12.2 Recommended accessories

(not included in delivery)

Part	Order number
DrillCut-X® ARTHRO shaver handpiece up to 8000 rpm, for use with UNIDRIVE® S III ARTHRO SCB, as of software status 1.10. The handpiece must NOT be treated with universal spray!	28200 DX
Valve body, with valve lock and simple lever, for use with DrillCut-X® ARTHRO shaver handpiece 28200 DX	28200 DC
Valve body, without valve lock and simple lever, for use with DrillCut-X® ARTHRO shaver handpiece 28200 DX	28200 DA
Valve lock, for valve body 28200 DC, for use with DrillCut-X® ARTHRO shaver handpiece 28200 DX	28200 DB
Shaver handpiece 8K	287210 35
Shaver handpiece 8K	287210 45
Shaver handpiece 15K	287210 30
Hybrid shaver handpiece	287210 50
Handpiece irrigation adaptor	28205 SA
Outer blade irrigation adaptor	28205 SAA
Inner blade irrigation adaptor	28205 SAI
Small joint handpiece 6K	287210 40
Dual-pedal footswitch	200168 31
Triple-pedal footswitch	200128 32
Multifunction handpiece	287210 36
Jacobs chuck 4 mm	287210 31
Jacobs chuck 6.5 mm	287210 37
Jacobs chuck (without key) 6 mm	287210 38
Sagittal Saw	287210 39
Pin driver	287210 33
Synthes style driver	287210 34

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**12.3 Recommended shaver blades
and burrs**

Aggressive cutter,
working length 180 mm, with keyhole, for removal
of degenerative meniscus tissue and adhesion
filaments.

Item	Cat. no.
------	----------

*Sterile, single use,
pack of 6*

Diameter 4.5 mm (violet), color code green	28208 BKS
---	-----------

Full radius resector,

working length 180 mm, with keyhole, for partial,
subtotal and total synovectomy, for removal of
degenerative meniscus tissue and adhesion fila-
ments.

Item	Cat. no.
------	----------

*Sterile, single use,
pack of 6*

Diameter 4.5 mm (violet) color code yellow	28208 CCS
Diameter 5.5 mm (dark gray), color code yellow	28208 CDS

Aggressive full radius resector,

working length 180 mm, with keyhole, for par-
tial, subtotal and total synovectomy, for removal
of degenerative meniscus tissue and adhesion
filaments.

Item	Cat. no.
------	----------

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28208 EGS

Curved full radius resector,
working length 180 mm, with keyhole, distal end
angled at 15°, for partial, subtotal, and total syno-
vectomy, for removal of degenerative meniscus
tissue and adhesion filaments. The curvature
enables inaccessible meniscus and cartilage
areas, and also synovial areas to be reached.

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28208 EHS

Curved aggressive full radius resector,
working length 180 mm, with keyhole, distal end
angled at 15°, for partial, subtotal, and total syno-
vectomy, for removal of degenerative meniscus
tissue and adhesion filaments. The curvature
enables inaccessible meniscus and cartilage
areas, and also synovial areas to be reached.

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28208 IDS

Semi-hooded barrel burr,
working length 180 mm, for removing large areas,
for removing osteophytes or a notchplasty; for
smoothing bone irregularities, for milling the sub-
chondral bone during subchondral abrasion for
induction of fibrous cartilage, extensive cartilage
ablation.

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Pro Line Shaver Blades,
working length 120 mm. Partial, subtotal and total synovectomy, removal of degenerative meniscus tissue, removal of cicatricial strands and adhesion filaments. Removal of osteophytes or a notchplas- ty; smoothing of bone irregularities, milling of the subchondral bone during subchondral abrasion for induction of fibrous cartilage, removal of large areas of cartilage, e.g., with subacromial decom- pression, for safe removal of meniscus tissue and adhesion filaments.

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Aggressive Pro Line Shaver Blades,

working length 120 mm. Partial, subtotal and total synovectomy, removal of degenerative meniscus tissue, removal of cicatricial strands and adhesion filaments. Removal of osteophytes or a notchplasty; smoothing of bone irregularities, milling of the subchondral bone during subchondral abrasion for induction of fibrous cartilage, removal of large areas of cartilage, e.g., with subacromial decompression, for the safe removal of meniscus tissue and adhesion filaments.

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Aggressive cutter,
working length 120 mm, for removal of degenerative meniscus tissue and filaments of adhesion.

Item	Cat. no.
------	----------

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End cutter,
working length 120 mm, for removal of degenerative meniscus tissue, cicatrized areas of synovial membrane; for removing large flaps of cartilage.

Item	Cat. no.
------	----------

Reusable, sterilizable

Diameter 3.5 mm, color code red	28205 BB
Diameter 4.2 mm, color code red	28205 BC
Diameter 5.5 mm, color code red	28205 BD

Sterile, single use, pack of 6

Diameter 3.5 mm, color code red/light gray	28205 BBS
Diameter 4.2 mm, color code red/middle gray	28205 BCS
Diameter 5.5 mm, color code red/dark gray	28205 BDS

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Full radius resector,
working length 120 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
------	----------

Reusable, sterilizable

Diameter 3.5 mm, color code yellow	28205 CB
Diameter 4.2 mm, color code yellow	28205 CC
Diameter 5.5 mm, color code yellow	28205 CD

*Sterile, single use,
pack of 6*

Diameter 3.5 mm, color code yellow/light gray	28205 CBS
Diameter 4.2 mm, color code yellow/middle gray	28205 CCS
Diameter 5.5 mm, color code yellow/dark gray	28205 CDS

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Aggressive full radius resector,
working length 120 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
------	----------

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Round burr,
working length 120 mm, for removal of osteo-
phytes, for smoothing bone irregularities, for
milling the subchondral bone during subchondral
abrasion, for induction of fibrous cartilage.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 4.2 mm, color code orange	28205 FC
Diameter 5.5 mm, color code orange	28205 FD
Diameter 6.5 mm, color code orange	28205 FE
<i>Sterile, single use, pack of 6</i>	
Diameter 4.2 mm, color code orange/middle gray	28205 FCS
Diameter 5.5 mm, color code orange/dark gray	28205 FDS
Diameter 6.5 mm, color code orange/black	28205 FES

**Finish barrel burr
(10 flutes, 4 notches),**
working length 120 mm, for removing large areas,
e.g. during subacromial decompression at the
shoulder; for removing osteophytes or a notch-
plasty; for smoothing bone irregularities, for milling
the subchondral bone in subchondral abrasion, for
induction of fibrous cartilage; for removing large
areas of cartilage, e.g. in subacromial decompres-
sion.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 4.2 mm, color code purple	28205 GD
Diameter 5.5 mm, color code purple	28205 GE

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**Aggressive barrel burr
(8 flutes, 4 notches),**
working length 120 mm, for removing large areas, e.g. during subacromial decompression at the shoulder; for removing osteophytes or a notch-plasty; for smoothing bone irregularities, for milling the subchondral bone in subchondral abrasion, for induction of fibrous cartilage; for removing large areas of cartilage, e.g. in subacromial decompression.

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Slotted Whisker,
working length 120 mm, for safe removal of
meniscus tissue and filaments of adhesion.

Item Cat. no.

Reusable, sterilizable

Diameter 4.2 mm,
color code light orange 28205 LC

Diameter 5.5 mm,
color code light orange 28205 LD

*Sterile, single use,
pack of 6*

Diameter 4.2 mm,
color code light orange/middle gray 28205 LCS

Diameter 5.5 mm,
color code light orange/dark gray 28205 LDS

Auger,
working length 120 mm, for removing large areas,
e.g., during subacromial decompression at the
shoulder; for removing osteophytes or a notch-
plasty; for smoothing bone irregularities, for milling
the subchondral bone in subchondral abrasion to
induce fibrous cartilage; for removing large areas
of cartilage, e.g. in subacromial decompression;
for safe removal of meniscus tissue and filaments
of adhesion.

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**12.4 Recommended shaver blades
and burrs for small joint
handpiece**

Aggressive cutter,

working length 70 mm, for removal of degenerative meniscus tissue and filaments of adhesion.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code green	28206 AA
Diameter 3.5 mm, color code green	28206 AB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code green/white	28206 AAS
Diameter 3.5 mm, color code green/light gray	28206 ABS

Full radius resector,

working length 70 mm, for partial, subtotal, and total synovectomy, for removal of degenerative meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code yellow	28206 CA
Diameter 3.5 mm, color code yellow	28206 CB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code yellow/white	28206 CAS
Diameter 3.5 mm, color code yellow/light gray	28206 CBS

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Aggressive full radius resector,
working length 70 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code blue	28206 DA
Diameter 3.5 mm, color code blue	28206 DB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code blue/white	28206 DAS
Diameter 3.5 mm, color code blue/light gray	28206 DBS

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Slotted Whisker,
working length 70 mm, for safe removal of menis-
cus tissue and filaments of adhesion.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code light orange	28206 LA
Diameter 3.5 mm, color code light orange	28206 LB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code light orange/white	28206 LAS
Diameter 3.5 mm, color code light orange/light gray	28206 LBS

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Small Joint Burr,
working length 70 mm, for removal of osteo-
phytes, smoothing of bone irregularities.

Item	Cat. no.
------	----------

*Sterile, single use,
pack of 6*

Diameter 2.5 mm (white), color code orange	28206 FAS
---	-----------

Diameter 3.5 mm (light gray), color code orange	28206 FBS
--	-----------

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**12.5 Accessories for shaver
handpieces and SL blades**

Irrigation adaptor, for validated cleaning of
shaver handpieces 28205 SA

Irrigation adaptor, for validated cleaning of
SL outer blades 28205 SAA

Irrigation adaptor, for validated cleaning of
SL inner blades 28205 SAI

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12.6 Accessories for pin driver

Drill wire, diameter 2.4 mm,
spiral-ground 28729 D

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Desgleichen, Blattbreite 8 mm 28207 SLE

12.7 Accessories for sagittal saw

Saw blade, blade width 10 mm,
blade thickness 0.4 mm, cutting depth 0.6 mm,
working length 50 mm 28207 SLA

Same, blade width 14 mm 28207 SLB

Saw blade, blade width 10 mm,
blade thickness 0.4 mm, cutting depth 0.6 mm,
with stopper at 15 mm,
working length 50 mm 28207 SLC

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Saw blade, blade width 12 mm,
blade thickness 0.4 mm, cutting depth 0.6 mm,
working length 22 mm 28207 SLF

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Desgleichen, Blattbreite 10 mm 28207 SLH

Saw blade, blade width 6 mm,
blade thickness 0.4 mm, cutting depth 0.6 mm,
working length 31 mm 28207 SLG

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13 Single-use blades and burrs

13.1 Intended use

The single-use sterile blades/burrs are intended for use with the UNIDRIVE® S III ARTHRO.

Use of the single-use sterile blades/burrs in any way other than that described is not permitted.

13.2 Use of the single-use products and explanation of the symbols employed

WARNING: Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the content of the pack must not be used under any circumstances.

WARNING: This product is intended for **single use** and must not be reused. Discard after use.
Do not allow the blade to run 'dry'. The recommended speed of rotation for the blades is 4,000 rpm.

NOTE: After use, dispose of the product according to instructions. National laws and regulations must be observed.

NOTE: This product is individually wrapped in protective packaging and sterilized with gamma radiation.

Manufacturer:

Manufacturer in accordance with Directive 93/42/EEC:

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Traceability

To ensure traceability, the self-adhesive identification labels enclosed in the original packaging can be kept and stuck in the operation report and in the patient dossier.

The person responsible for stores at the hospital can maintain a record of the products currently in store, including article no. (REF) and lot no. (LOT) so that rapid action can be taken in the event of changes or recalls.

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Electromagnetic Compatibility (EMC) Information

14 Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this section during installation and commissioning.

The UNIDRIVE® S III ARTHRO model 28 7230 20-1 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The UNIDRIVE® S III ARTHRO model 28 7230 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this section provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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**Electromagnetic
Compatibility (EMC)
Information**

WARNING: The UNIDRIVE® S III ARTHRO model 28 7230 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the UNIDRIVE® S III ARTHRO model 28 7230 20-1.

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Table 1		
Guidance and manufacturer's declaration – electromagnetic emissions		
UNIDRIVE® S III ARTHRO model 28 7230 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III ARTHRO model 28 7230 20-1 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group1	UNIDRIVE® S III ARTHRO model 28 7230 20-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The UNIDRIVE® S III ARTHRO model 28 7230 20-1 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions acc. to CISPR 11	Class B	
Harmonic emissions acc. to IEC 61000-3-2	Class D (see NOTE below)	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	not applied (see NOTE below)	
 NOTE: The measurement of harmonic emissions is not required (power < 75 W) The measurement of voltage fluctuations/flicker is not required (activated manually, unlikely to produce significant voltage fluctuations or flicker)		

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
UNIDRIVE® S III ARTHRO model 28 7230 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III ARTHRO model 28 7230 20-1 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	1 kV/2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	0 V, 10 msec 92 V, 100 msec 168 V, 500 msec 0 V, 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
UNIDRIVE® S III ARTHRO model 28 7230 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III ARTHRO model 28 7230 20-1 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF disturbance acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _m	Portable and mobile RF communications equipment should be used no closer to any part of the UNIDRIVE® S III ARTHRO model 28 7230 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF disturbance acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 6 Recommended separation distances between portable and mobile RF communications equipment and the UNIDRIVE® S III ARTHRO model 28723020-1			
UNIDRIVE® S III ARTHRO model 28723020-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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KARL STORZ Endoscopy Canada Ltd.
1711 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseal.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseal.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326 - Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4, kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TD, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch.
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch.
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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**Important information for
users of KARL STORZ
devices and instruments**

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**1 Important information for
users of KARL STORZ
devices and instruments**

**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

This instruction manual is designed to provide help with the proper use, cleaning, and disinfection of the POINT SETTER®. For this reason, please read these instructions carefully and keep them as a ready reference.

MITAKA KOHKI and KARL STORZ are constantly working on the further development of all products. Please understand that changes to the scope of supply, design, equipment, and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

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
- ① Safety lock knob
- ② Pneumatic joints/swivel axis
- ③ Release switch
- ④ Manual ball joint
- ⑤ Mount for clamping jaws
- ⑥ Mount for holding arm

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
- ⑦ Swivel joint
(swivelling the arm around a joint)
- ⑧ Lock
(moving the cylinder rail to right/left)
- ⑨ Swivel joint
(turning the cylinder rail around an axis)
- ⑩ Lock
(fixing to operating table rail)

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1	Important information for users of KARL STORZ devices and instruments.....	III
2	Image of the unit.....	IV
2.1	MITAKA KOHKI POINT SETTER®	IV
2.2	Holding arm	V
2.3	Operating table adaptor	VI
3	Safety instructions.....	4
3.1	Explanation of warnings and cautions.....	4
3.2	General safety information	5
3.3	Intended use.....	7
3.3.1	Safety device	7
3.3.2	Contraindications.....	7
3.3.3	Safety precautions for the use of the MITAKA KOHKI POINT SETTER®	8
3.3.4	User qualification	8
4	Description of the device.....	9
4.1	Storage conditions.....	9
4.2	Warning signs on the unit	10
4.3	Explanation of symbols	10
5	Initial use	11
5.1	Unpacking	11
5.2	General information on assembly.....	11
6	Operating table adaptor.....	14
6.1	Fastening the adaptor to the operating table.....	14
6.2	Positioning the operating table adaptor	14
6.2.1	Turning the cylinder rail	15
6.2.2	Moving the cylinder rail	16
6.2.3	Turning the long swivel arm	17
6.2.4	Turning the short swivel arm	17

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7	Holding arm	18
7.1	Fastening the holding arm to the operating table adaptor	18
7.2	Connecting up the compressed air hose of the holding arm	18
7.3	Positioning the holding arm	20
7.4	Sterile covering of the MITAKA KOHKI POINT SETTER®	22
7.5	Fixing instruments and telescopes to the holding arm	22
7.5.1	Use of KARL STORZ clamping jaws with the KSLOCK adaptor 28272 UA	22
7.5.2	Use of the MITAKA KOHKI clamping jaws.....	24
7.6	Positioning the manual ball joint	25
7.7	Functional test before surgery.....	26
8	Reprocessing	27
8.1	Wipe-down disinfection of the MITAKA KOHKI POINT SETTER®	28
8.1.1	Removing the clamping jaw from the holding arm	28
8.1.2	Performing a wipe-down disinfection	28
8.2	Cleaning, disinfection, and sterilization of the clamping jaws	29
8.2.1	Manual cleaning.....	30
8.2.2	Manual disinfection	31
8.2.3	Assembly, testing, care and packaging	31
8.2.4	Steam sterilization.....	32
8.2.5	Pre-vacuum method	32
8.2.6	Reprocessing overview table	34
9	Maintenance	36
9.1	Holding arm.....	37
9.2	Operating table adaptor.....	39

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3 Safety instructions

3.1 Explanation of warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this instrument. To make the signal words stand out more clearly, they are accompanied by a pictogram.

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CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the instrument.

NOTE: A Note indicates special information about the operation of the instrument, or clarifies important information.

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3.2 General safety information

 **WARNING:** Risk of injury/material damage if the instruction manual/warnings are not observed: Read and observe the instruction manual and warnings for the devices and the devices used in combination.

 **WARNING:** Safety instruction: This product must only be operated by a trained member of medical personnel who has been correctly instructed in its use.

 **WARNING:** Risk of injury to the patient/user: Damaged or incomplete units/accessories can cause malfunctions. Always check units/accessories before and after use to ensure that nothing is missing/damaged. Replace incomplete/damaged units/accessories and do not continue to use them.

 **WARNING:** Risk of injury to the patient/user: Malfunctions can result in improper functioning. Perform a test for proper functioning before every use.

 **WARNING:** Repair, conversion, maintenance, and other technical services must be performed exclusively by an appointed service person of MITAKA KOHKI or technicians trained by MITAKA KOHKI. MITAKA KOHKI's warranty will not cover any mechanical failure or liability caused by failure to comply with this instruction manual.

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WARNING: Risk of injury/material damage: All parts requiring the attention of technical services or repair must be replaced with MITAKA KOHKI products only.

WARNING: In the USA, Federal law restricts this device to sale or use by or on the order of a physician.

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WARNING: Risk of fatal injury to the patient if the working load of 2 kg is exceeded.

WARNING: Risk for patient: The POINT SETTER® should be serviced every month as described in the Maintenance chapter of this manual.

WARNING: Risk for patient: If operating pressure is too low (< 6 bar), the arm will be stiff and hard to move.

WARNING: Risk for patient: If operating pressure is too high (> 8 bar), the arm could spontaneously unlock or refuse to lock when the switch is released.

CAUTION: Risk of material damage: Tighten locks and joints manually and not with tools, as they may otherwise be damaged.

NOTE: The POINT SETTER® was designed according to standard IEC 60601-1/EN 60601-1.

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3.3 Intended use

The POINT SETTER® is a mechanical pneumatic holding arm which is affixed to an operating table using an operating table adaptor. The holding arm is used to position and hold optical systems (endoscopes, VITOM® incl. camera), retractors, or instruments.

3.3.1 Safety device

The POINT SETTER® has a reliable safety mechanism, which locks the holding arm if the power fails or gas pressure is too low.

3.3.2 Contraindications

No contraindications directly related to the product are currently known. Use of the unit is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through its use, e.g., due to the general condition of the patient, or if the endoscopic method itself is contraindicated.

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3.3.3 Safety precautions for the use of the MITAKA KOHKI POINT SETTER®

The POINT SETTER® must be used according to the medical rules and procedures of endoscopy recognized for this method.

3.3.4 User qualification

The POINT SETTER® may only be used by persons with an appropriate medical qualification and who are acquainted with the requisite technique.

The information given in this manual only serves to instruct in the correct handling, cleaning, and disinfection of the POINT SETTER® and the sterilization of the accessories. It is not intended as an introduction to the technique.

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4 Description of the device

The POINT SETTER® is intended as a holding arm for optical systems (endoscopes, VITOM® incl. camera), retractors, or instruments, which can be positioned and then held in position. The POINT SETTER® consists of a holding arm and an operating table adaptor. The holding arm can be locked and released using a release switch. The holding arm can be attached on either side rail of the operating table using an operating table adaptor.

The operating table adaptor can be adjusted to allow the holding arm to be at the best position for the surgeon. The holding arm has a reliable safety mechanism, which locks the arm if the power fails or gas pressure is too low.

4.1 Storage conditions

CAUTION: Risk of material damage: Observe the following storage conditions for the unit and the accessories.

- Store the unit in a dry place.
- Protect the unit from water splashes and other liquids.
- Do not store the unit in a place susceptible to vibration or shock.
- The storage temperature should be room temperature (approx. 20 °C).
- Store the cables and hoses together after cleaning.

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


Hersteller

4.2 Warning signs on the unit

Should any of the warning signs be missing from the positions specified, contact the distributor to rectify the problem. Do not use the unit during this time.

POINT SETTER®
MITAKA KOHKI CO., LTD
1-18-8, Nozaki, MITAKA-SHI,
TOKYO 181-0014 JAPAN

Only to be operated by trained personnel!

Using gas pressure 6.0-8.0 kgf/cm²

Max. load 2 kg.

Ensure that all parts are securely fastened.

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5 Initial use

WARNING: Risk of infection: The operating table adaptor and the holding arm are delivered unsterilized. Wipe-down disinfect the operating table adaptor and the holding arm before the first use and before every subsequent use to prevent infections.

WARNING: Risk of infection: The accessories, e.g., clamping jaws, are delivered unsterilized. Clean, disinfect, and sterilize the clamping jaws before initial use and before every subsequent use using validated procedures to prevent infection.

5.1 Unpacking

Carefully unpack the POINT SETTER® and remove the unit and its accessories from their packaging. Check the delivery for missing items and possible damage.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

5.2 General information on assembly

WARNING: Risk of injury/Risk of infection: Repositioning the holding arm during the surgical procedure may be difficult and it may prove necessary to have an unsterile person to manipulate the operating table adaptor. It is thus important to test the operating range of the POINT SETTER® with equipment before commencing surgery and make any necessary changes in advance.

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WARNING: Risk of injury to the patient due to unstably fixed POINT SETTER®. Some operating tables are curved by design. The design features of the respective operating table must be taken into account when operating the POINT SETTER®.

WARNING: Risk of injury to the patient due to unstable POINT SETTER®. The stability of the POINT SETTER® is dependent upon the stability of the operating table's side rails. Check the operating table thoroughly. You must not use the POINT SETTER® if anything is loose on the operating table, e.g., screws.

The POINT SETTER® holding arm is fixed to the side rails of the operating table using the operating table adaptor.

The POINT SETTER® holding arm and the operating table adaptor can be very flexibly adjusted, which enables the exact positioning of the instrument to be held.

The position of the POINT SETTER® arm on the operating table should be selected carefully in accordance with the surgical intervention to be carried out.

Make the following preparations in particular:

- The POINT SETTER® holding arm should be attached before covering the patient with sterile drapes.
- The POINT SETTER® holding arm should be positioned to suit the intervention and the position of the patient. The position should be selected in such a way that the operating field can be reached comfortably with the arm.

WARNING: Danger to the patient: The position at which the POINT SETTER® arm is fixed to the operating table should not be changed after sterile draping and commencement of the surgery.

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
- Locate the holding arm in the best position for the surgical procedure, allowing for use of C-arms or other surgical instruments.

The individual steps are explained below.

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- ① NOTE:** The brake knobs on the joints are adjustable and, if necessary, can be turned to another position without opening the lock or the swivel joint.

Moving the brake knob

1. Pull the knob of the brake knob towards you.
2. Turn the knob to the required position.
3. Release the knob.

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

4

6 Operating table adaptor

6.1 Fastening the adaptor to the operating table

1. Engage the adaptor on the side rail of the operating table.
2. Firmly tighten the operating table rail lock.

6.2 Positioning the operating table adaptor

You can use different adjustment mechanisms (manual) to position the operating table adaptor as required for the particular surgical procedure.

WARNING: Risk of fatal injury to the patient from inadvertent movement of the POINT SETTER® as a result of excessive weight: The weight of the holding arm prohibits the safe adjustment of the operating table adaptor. When the holding arm is attached, do not attempt to adjust the operating table adaptor.

WARNING: Risk of fatal injury to the patient if the operating table adaptor is inadvertently moved. Ensure that the toothing on the swivel joints interlocks correctly when closing, and that all joints and locks are closed tightly.

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
6.2.1 Turning the cylinder rail

1. Hold the operating table adaptor in position.
2. Open the swivel joint by turning the brake knob counterclockwise (Image 5).

3. Turn the cylinder rail to the required position.

! WARNING: Risk of injury from incorrectly closed locks and joints: Ensure that the toothing on the swivel joints interlocks correctly when closing, and that all joints and locks are closed tightly.

4. Close the swivel joint by turning the brake knob clockwise.

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
9

6.2.2 Moving the cylinder rail

WARNING: Risk of injury from incorrectly closed locks and joints: Ensure that the toothing on the swivel joints interlocks correctly when closing, and that all joints and locks are closed tightly.


1. Unscrew and remove the end cap (Image 8).
2. Open the cylinder rail lock by turning the respective brake knob counterclockwise.

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
NOTE: If necessary, the cylinder rail can also be pushed in from the opposite side.

3. Push the cylinder rail to the required position.

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4. Close the cylinder rail lock by turning the respective brake knob clockwise (Image 11).

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6.2.3 Turning the long swivel arm

1. Open the swivel joint by turning the thumb screw counterclockwise.
2. Turn the long swivel arm to the required position.

WARNING: Risk of injury from incorrectly closed locks and joints: Ensure that the toothing on the swivel joints interlocks correctly when closing, and that all joints and locks are closed tightly.

3. Close the swivel joint by turning the thumb screw clockwise.

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6.2.4 Turning the short swivel arm

1. Open the swivel joint by turning the thumb screw counterclockwise.
2. Turn the short swivel arm to the required position.

WARNING: Risk of injury from incorrectly closed locks and joints: Ensure that the toothing on the swivel joints interlocks correctly when closing, and that all joints and locks are closed tightly.

3. Close the swivel joint by turning the thumb screw clockwise.

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
7 Holding arm

7.1 Fastening the holding arm to the operating table adaptor

1. Open the mount for the holding arm by turning the thumb screw counterclockwise.
2. Insert the foot of the holding arm into the groove (Image 15).
3. Fix the holding arm in the mount by turning the thumb screw clockwise.
4. Check the positioning and retighten the thumb screw if necessary.

7.2 Connecting up the compressed air hose of the holding arm

WARNING: Risk of injury: Positioning via the pneumatic joints of the holding arm only functions when the compressed air supply is on.

WARNING: Risk of injury: Excessive pressure is extremely dangerous. Nitrogen gas or compressed air pressure must not exceed 6-8 bar or 90-115 psi.

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WARNING: Risk for patient if the performance of the holding arm is impaired when the pressure is too low: If the pressure is too low (<6 bar; 90 psi), the joints of the holding arm do not open completely. The holding arm is difficult to move in this case. Before the surgery, check the holding arm for proper function and mobility, and ensure that the operating pressure is 6-8 bar or 90-115 psi.

1. Ensure that the operating pressure is 6-8 bar or 90-115 psi.
2. Connect the hose of the holding arm to the hospital's compressed air supply (Image 16).

NOTE: If the central compressed air supply delivers a pressure > 8 bar, the pressure limiter 2761291 can be used. The pressure limiter is inserted between the connecting tube for the compressed air supply and the POINT SETTER® arm.

NOTE: When using the compressed air limiter 2761291, follow the accompanying instruction leaflet 96216314DF.

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
7.3 Positioning the holding arm

You can use different adjustment mechanisms (manual and pneumatic) to position the holding arm as required for the particular surgical procedure.

WARNING: Risk of fatal injury to the patient as a result of unstable operating table adaptor or incorrect mounting:

As the holding arm is attached to the operating table adaptor and the latter is fastened to the side rails of the operating table, the stability of the holding arm is directly dependent on the stability of the side rails of the operating table and the correct mounting of the operating table adaptor. Ensure that the operating table, the side rails and the operating table adaptor are stable and not damaged. Check the correct mounting of all components.

WARNING: Risk of fatal injury to the patient if the operating table adaptor is inadvertently moved. Ensure that the toothings on the swivel joints interlocks correctly when closing, and that all joints and locks are closed tightly.

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
WARNING: Risk of fatal injury to the patient due to inadvertent activation of the release switch and movement of the holding arm. After positioning the holding arm, it is essential to set the safety lock knob on the release switch to the 'Lock' position and check its positioning.

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CAUTION: Risk of material damage: Positioning via the pneumatic ball joints or the swivel joint only functions when the compressed air supply is switched on.

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1. Set the safety lock knob on the release switch to the 'Free' position.
2. Hold the holding arm tightly.
3. Press the release switch and hold it depressed (Image 18).
4. Place the holding arm in the required position.
5. Release the release switch.
6. Set the safety lock knob back to the 'Lock' position (Image 19).
7. Check the correct positioning of the holding arm.

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7.4 Sterile covering of the MITAKA KOHKI POINT SETTER®

Putting on the sterile cover

WARNING: Risk of infection: The POINT SETTER® is not sterile. Cover the POINT SETTER® carefully with a sterile cover.

WARNING: Risk for the patient if the holding arm's range of movement is hindered because the sterile cover is too tight: Do not apply the sterile cover too tightly. Check the range of movement before fastening the sterile cover.

1. Starting at the distal end of the holding arm, pull the sterile cover carefully over the entire POINT SETTER® (Image 20).

Fastening the sterile cover

1. Fasten the sterile cover using the adhesive tapes attached at various points (Image 21).

7.5 Fixing instruments and telescopes to the holding arm

7.5.1 Use of KARL STORZ clamping jaws with the KSLOCK adaptor 28272 UA

As described, cover the POINT SETTER® holding arm with a sterile cover.

1. Insert the adaptor 28272 UA into the mount provided for this purpose on the end of the holding arm (Image 22).

2. Tighten up the screw for the clamping jaw mount (Image 23).

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
Fixing the clamping jaw to the KSLOCK quick-release coupling

1. Release the thumb screw of the quick-release coupling (female) on the KSLOCK adaptor by turning it as far as it will go ("OPEN" position).
2. Select the appropriate clamping jaw for the instrument.
3. Attach the clamping jaw to the KSLOCK adaptor so that it clicks into place. Secure the quick-release coupling by turning it as far as it will go ("CLOSE" position).

NOTE: When locked correctly, the pin of the quick-release coupling is not visible (Image 25). If the pin of the quick-release coupling is visible (Image 26), please repeat step 3.

Inserting the instrument/telescope into the clamping jaw

NOTE: When fastening rigid telescopes to the optical block, the correct positioning of the clamping jaw must be ensured. With large telescope heads (10 mm telescopes), the clamping jaw 28272 UGN must be mounted flush with the upper edge of the telescope head (Image 27).

NOTE: When using rigid 10 mm telescopes, the clamping cylinder 28272 C or 28272 CN can also be mounted on the shaft of the telescope. The telescope with the clamping cylinder is then held by the clamping jaw 28272 UGK.

WARNING: Risk of fatal injury to the patient if the total load of 2 kg is exceeded. When new, the POINT SETTER® can hold telescopes, endocameras and instruments up to a total load of 2 kg. Do not exceed the maximum permissible load.

WARNING: Risk of fatal injury to the patient: Over time the brakes will weaken. Always test the arm with the load presented by the instrument it is intended to hold before using it on a patient. Test the function of the brakes as described in chapter 9 of this manual.

1. Insert the instrument or telescope into the clamping jaw, put it in the required position and fix it there using the wheel.

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7.5.2 Use of the MITAKA KOHKI clamping jaws

NOTE: The clamping jaw goes over the sterile cover. This clamps the sterile cover.

WARNING: Risk of infection from damaged sterile cover: Do not pull on the sterile cover when fastening the clamping jaw as this can cause damage.

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1. Insert the clamping jaw into the mount provided for this purpose on the end of the holding arm (Image 31).
2. Tighten up the screw for the clamping jaw mount (Image 32).

Fixing the endoscope to the mounted clamping jaw

1. Insert the sterile endoscope into the sterile clamping jaw.
2. Tighten the screw (clamping jaw).

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7.6 Positioning the manual ball joint

NOTE: You can position the endoscope and the clamping jaw by adjusting the manual ball joint.

WARNING: Risk of fatal injury to the patient if the manual ball joint of the holding arm is inadvertently moved. Tighten the thumb screw firmly (manual ball joint) and check its correct positioning.

1. Open the manual ball joint by turning the thumb screw counterclockwise.
2. Position the endoscope with clamping jaw optimally.
3. Close the manual ball joint by turning the thumb screw clockwise.

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38

7.7 Functional test before surgery

WARNING: Risk of injury: Before using the POINT SETTER® on the patient, test the holding arm with the planned instrument (load) in order to get accustomed to the weight.

WARNING: Risk of injury: Activating the release switch releases all the pneumatic joints of the POINT SETTER® and the weight of the arm and the working load are placed in the hands of the user.

It is thus important to hold the handle tightly and prepare to take the weight to avoid uncontrolled dropping of the arm.

Once you have brought the arm into position, ensure that all pneumatic joints are completely locked before releasing the handle. Depending on the weight of the attached load, the arm may drop slightly.

WARNING: Risk of fatal injury to the patient due to inadvertent activation of the release switch and movement of the holding arm during the surgical procedure. After positioning the holding arm, it is essential to set the safety lock knob on the release switch to the 'Lock' position and check its positioning.

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8 Reprocessing

 **NOTE:** The 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions must be followed. The cleaning, disinfection, care and sterilization procedures are explained therein in detail.

 **NOTE:** An overview table detailing suitable reprocessing measures for individual items can be found on page 34.

 **WARNING:** Risk of infection: If these inspections and cleaning regulations are not observed, there is a risk of contamination.

 **CAUTION:** Risk of material damage: Liquid must not be allowed to enter the holding arm's joints under any circumstances.

 **NOTE:** Throughout this sterilization section, the reusable sterilizable accessories that are attached to the holding arm outside the sterile cover, are collectively referred to as 'clamping jaw' for purposes of simplification.

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8.1 Wipe-down disinfection of the MITAKA KOHKI POINT SETTER®

WARNING: Wipe-down disinfect the operating table adaptor and the holding arm before the first use and before every subsequent use to prevent contact contamination.

8.1.1 Removing the clamping jaw from the holding arm

1. Undo the screw.
2. Remove the holder with the attached endoscope from the mount.

8.1.2 Performing a wipe-down disinfection

CAUTION: Risk of material damage: Do not clean the POINT SETTER® in a washer/disinfector and do not sterilize the POINT SETTER®.

1. Perform a wipe-down disinfection (no solutions predominantly containing alcohol) of the exterior surfaces of the POINT SETTER®.
We exclusively recommend using products based on ammonium chloride or quaternary ammonium chloride.

NOTE: Always observe the chemical manufacturers' specifications regarding application and concentration.

WARNING: Risk of injury: The pneumatic joints of the holding arm must not be oiled or greased under any circumstances, as this reduces the holding force. We recommend degreasing the pneumatic joints with alcohol or acetone, if necessary, at regular intervals.

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8.2 Cleaning, disinfection, and sterilization of the clamping jaws

 **WARNING:** Risk of infection: The accessories, e.g., clamping jaws, are delivered unsterilized. Clean, disinfect, and sterilize the clamping jaws before initial use and before every subsequent use using validated procedures to prevent infection.

 **WARNING:** Risk of infection/material damage: Observe the following instructions carefully when cleaning, disinfecting, and sterilizing the clamping jaw.

 **WARNING:** Risk of infection/material damage: Any deviations from the recommended parameters for cleaning and sterilization should be validated by the user.

 **CAUTION:** Risk of material damage: The POINT SETTER® arm, operating table adaptor, connecting tubes and pressure limiter cannot be sterilized. In these cases, use a wipe-down disinfection only.

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8.2.1 Manual cleaning

CAUTION: Risk of infection: Wear protective gloves, clothing, and a face mask when cleaning a contaminated clamping jaw.

1. Thoroughly rinse heavy soiling off all parts with completely demineralized water.

NOTE: Manual cleaning is recommended. However, areas which are difficult to reach can be cleaned using an ultrasound bath. It is important to follow the manufacturer's instructions for operating the ultrasound bath.

2. Immerse the clamping jaw fully in cleaning solution. An enzymatic or pH-neutral cleaning solution is recommended. The use of an alkaline cleaning agent can lead to cosmetic changes. However, these changes do not have a detrimental effect on functioning and safety.
3. Remove any residual blood, protein material and contaminants from the clamping jaw with a disinfected brush or sponge, soft cloths or a cotton tip applicator.
4. Clean inaccessible parts with disinfected cleaning brushes of a suitable size.
5. Triple rinse the clamping jaw with completely demineralized water for a minimum of one minute for each rinse.
6. The rinse water should be discarded at the end of each rinse, as it will then be contaminated.

NOTE: National laws and regulations must be observed.

7. Dry the clamping jaw with a low-lint soft cloth or sterile compressed air.

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8.2.2 Manual disinfection

The instrument must be opened/dismantled as far as possible and then fully immersed in the disinfectant solution. Ensure that all external/internal surfaces and surfaces with restricted access are completely soaked. Following the exposure time, the instrument must be rinsed off/through several times with completely demineralized water or microbiologically pure/sterile water in order to remove all chemical residues. All of the surfaces, joints, openings, channels, and lumina are then dried completely, preferably with sterile compressed air. The compressed air gun with accessories (Art. No. 27660) is ideal for this purpose.

8.2.3 Assembly, testing, care and packaging

The cleaned/disinfected instrument must be visually inspected to ensure that it is clean, complete, dry and has not been damaged. If residues or contamination are still present, then the instrument must be manually cleaned and subjected to a full cleaning and disinfection procedure once more. A damaged or corroded instrument must be disposed of. A dismantled instrument must be assembled. Joints, threads and sliding surfaces must be lubricated locally with instrument oil/special grease.

Afterwards, carry out a functional check.

- ① **NOTE:** The oil used for this purpose must be suitable for the subsequent sterilization process (silicone-free and paraffin- or white oil-based).
- ① **NOTE:** During care procedures, use items from the catalog "Care, Sterilization and Storage Techniques".
- ① **NOTE:** Greased components (e.g., stop-cock) must be dismantled for sterilization so that steam penetration can be ensured.

Only standardized and approved packaging materials or systems must be used.

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8.2.4 Steam sterilization

WARNING: Risk of infection/material damage: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and calibrated.

WARNING: Risk of infection/material damage: Deviations from the recommended sterilization parameters must be validated by you as the user.

1. Place the clamping jaw in a sterilization tray.

MITAKA KOHKI has validated the following steam sterilization methods:

8.2.5 Pre-vacuum method

Pre-vacuum or high vacuum sterilization consists of four basic phases:

- conditioning phase;
- exposure phase;
- evacuation phase;
- drying phase.

The conditioning phase removes air from the chamber by creating a vacuum and then heats the clamping jaw by injecting steam. Sterilization occurs during the exposure phase, when the chamber reaches a temperature of 132.2 °C and pressure of 1.85 bar. The exposure phase in a pre-vacuum type of sterilizer last for 4.0 minutes. The evacuation phase removes the steam from the chamber. A new vacuum is created in the drying phase. The clamping jaw dries in approx. 20 minutes in a vacuum. The following conditions have been used to validate sterilization procedures in a pre-vacuum sterilizer:

Temperature:	132.2 to 133.3 °C
Pressure:	1.85 bar or 27 psi
Exposure time:	4 mins

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8.2.6 Reprocessing overview table

			Wipe-down disinfection	Preliminary cleaning				Cleaning				Sterilization											
				manual				by machine				STERRAD®											
Article	Description	Disassembly	Wipe-down disinfection	Immersion in cold water	Brushing the surfaces	Brushing the lumens	Rinsing out with a water gun	Ultrasonic treatment	Manual cleaning	Manual disinfection	Chemical disinfection	Thermal disinfection	Connector	Steam sterilization	STERRAD® 50, 100S, 200	STERRAD® NX	STERRAD® 100NX	STERRAD® 100NX 000 Cycle	AUSCUT® V-PRO™	Ethylene oxide (EO)	Steam & formaldehyde process, low temperature	STERIS® System 1	STERIS® System 1E
28272 UGN	Clamping jaw, large	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28272 UKN	Clamping jaw, small	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28272 UF	Clamping jaw, fiberscope	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28272 UA	POINT SETTER® KSLOCK adaptor	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28172 WR	Blade holder	-	x	x	x	x	x	-	x	x	-	-	-	-	-	-	-	-	-	-	-	-	-
28172 WX	Universal adaptor, 10-15 mm, pediatric	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28172 WV	Universal adaptor, 10-15 mm	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28172 WU	Universal adaptor, 5-10 mm	-	x	x	x	x	x	-	x	x	-	-	-	x	-	-	-	-	-	-	-	-	-
28172 P	POINT SETTER® arm	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28172 T	Table adaptor, simple	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28172 O	Table adaptor, original	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28272 WN	Connecting tube, Schrader	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28272 WA	Connecting tube, Drager	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28272 WB	Connecting tube, air-motor	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28272 WC	Connecting tube, compressor	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28272 WO	Connecting tube, open	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2761291	Pressure limiter	-	x	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

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9 Maintenance

 **NOTE:** The POINT SETTER® should be serviced monthly by trained technical staff. Test the POINT SETTER® as described in this chapter.

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9.1 Holding arm

1. Check the surfaces of the POINT SETTER®.

Expected result:

- No damage
- No contamination
- No loose parts

2. Check the compressed air hose and connections.

Expected result:

- Compressed air hoses are leak-proof
- Connections are leak-proof

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4. Check the brake holding force. The brakes on the POINT SETTER® must withstand a load of 2 kg and above; otherwise the unit should be sent to the factory for repair.

Expected result:

- The holding force of the brakes 2, 3, 4, and 5 is equivalent to 2 kg or more.

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10 Service and repair**10.1 Service**

For servicing, maintenance and any questions you may have in connection with the POINT SETTER® and accessories, contact the technical service.

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

10.2 Repair

Defective items of equipment may only be repaired by persons authorized by us and using original parts from MITAKA KONKI.

10.3 Repair program

During the repair period, you will generally be loaned a device, which must then be returned to KARL STORZ as soon as the repaired device is returned to you.

In Germany you can refer repairs to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries, please contact your local KARL STORZ subsidiary or authorized dealer.

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10.4 Important information

WARNING: Risk of infection from shipping of contaminated medical devices: In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

NOTE: We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes, or expansions which are not carried out by MITAKA KOHKI CO., LTD. or by experts authorized by MITAKA KOHKI CO., LTD. will invalidate all warranty rights.

MITAKA KOHKI CO., LTD. gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

10.5 Limitation of liability

As the manufacturer of this device, MITAKA KOHKI CO., LTD. shall only be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by MITAKA KOHKI CO., LTD.
- the electrical installation of the location in which the device is connected and operated conforms to the applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times
- the unit is tested as described in this manual.

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10.6 Warranty

The warranty period for the MITAKA KOHKI POINT SETTER® is 1 year. MITAKA KOHKI CO., LTD. guarantees the product against material and manufacturing defects. The warranty period, from the date of invoice, can be found in the respective instruction manual. Faulty material will be replaced free of charge during the warranty period provided that it has not been subject to improper use. This warranty does not cover wearing parts. In order not to void the warranty and to ensure that the unit functions properly, the information in the instruction manual must be followed and only accessories and spare parts from MITAKA KOHKI CO., LTD. may be installed and/or used.

 **NOTE:** The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):
KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

 **CAUTION:** Voiding of warranty: Opening of the unit or the performance of any repairs or modifications to the unit by unauthorized persons or modifications which do not comply with the applicable sections of the IEC 60601 set of standards and the EN ISO 10079-1 standard shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

Warranty rights which go beyond the described warranty, such as liability for consequential damage etc., have not been agreed.

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11 Technical description

For servicing, maintenance and any questions you may have in connection with the MITAKA KOHKI POINT SETTER® and accessories, contact the technical service.

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

11.1 Troubleshooting

Symptom:

The holding arm cannot be moved despite the release switch being pressed.

Remedy:

- Connect nitrogen hose.
- Use pressure of 6-8 bar.
- Turn the safety lock knob on the release switch to 'Free'.

Symptom:

Holding force too weak.

Remedy:

- Check weight of the instrument (maximum load is 2 kg).
- Check the joints of the holding arm.
- Check the operating table adaptor.

Symptom:

Endoscope or instrument slips.

Remedy:

- Check the individual clamps on the holding arm.
- Check the clamping jaw.
- Check the clamps on the operating table adaptor.
- Check the operating table.

NOTE: If the fault cannot be remedied, contact our local representative for service.

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11.4 Technical data

Holding arm	
Number of joints (incl. manual)	6
Arm length (nominal)	650 mm
Weight	6 kg
Maximum load	2 kg
Operating table adaptor	
Number of joints (slides)	3 (1)
Weight	6 kg

11.5 Storage and operating conditions

CAUTION: Risk of material damage: Observe the following storage conditions for the device and accessories.

- Store the device in a dry place.
- Protect the device from water splashes and other liquids.
- Do not store the device in a place susceptible to vibration or shock.
- Store the device at room temperature (approx. 20 °C).
- Store the cables and hoses together after cleaning.

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Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device. Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kazakhstan LLP
6, Saryarka str. BC "Arman", office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd
D-181, Okhla Industrial Area. Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

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KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company.
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016
People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company.
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

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
Important note for users of KARL STORZ devices and instruments

It is recommended to check the suitability of the products for the intended procedure prior to use.

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This instruction manual is intended to serve as an aid in the proper setup and installation of the KARL STORZ FI System. All of the necessary details and actions are clearly explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

General warning

Users of KARL STORZ devices should be thoroughly trained in their use and the applicable medical procedures. Use of these devices and accessories should not be undertaken until the user has familiarized himself/herself fully with the instruction manual, assembly, and care. Instruction manuals should be carefully studied and be available to the surgical team during the procedure. It is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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Fluorescence system video cart
(example)


- ① Video cart with isolation transformer
- ② KARL STORZ monitor
- ③ FI camera head or FI videoscope
- ④ PDD, AF or ICG telescope (not in figure)
- ⑤ FI camera control unit
- ⑥ D-LIGHT light source
- ⑦ Fluid light cables for rigid systems with proximal camera (not in figure)
- ⑧ Footswitch
- ⑨ SCB connecting cable

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**Important note for users of
KARL STORZ devices and
instruments**

General warningIII

**System images for
KARL STORZ
fluorescence systems**

Equipment cart (example)IV
TRICAM® system overview for
PDD, AF and ICG applications
with proximal camera heads3
TELECAM system overview for
PDD and AF applications
with videoscopes6
IMAGE 1 system overview for
ICG applications with proximal
camera head9

System description

Components of the KARL STORZ FI System 11
Combination of light cable
with endoscope 11

Safety instructions

Warnings and cautions 12
Note 12
Safety information for medical
electrical systems 13

Service and repair

Reprocessing 14
Preventative maintenance measures 14
Maintenance and safety check 15
Maintenance 15
Safety checks/
repeat inspections as per IEC 62353 15

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


TRICAM® system overview for PDD,
AF and ICG applications with
proximal camera heads

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


TELECAM system overview for PDD
and AF applications with
videoscopes

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

**IMAGE1 system overview for
ICG applications with proximal
camera head**

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Components of the KARL STORZ FI System

The KARL STORZ FI System comprises the KARL STORZ D-LIGHT C/AF/P system including suitable KARL STORZ fluid light cables, KARL STORZ FI endoscope or FI videoscope and KARL STORZ FI camera with coordinated KARL STORZ camera control unit.

All modifications or additions to this system require a complete new assessment and documentation for the purposes of IEC 60601-1-1 and IEC 60601-1.

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Combination of light cable with endoscope

For safety reasons, the diameters of fiberoptic/fluid light cables and endoscopes should be matched to each other.

NOTE: If there is an identification symbol on the light cable and the endoscope (see table), assign matching combinations according to the symbols shown.

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Overview: Recommended combinations*

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* Special endoscopes may deviate slightly

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Warnings and cautions

WARNING: All instruction manuals included with the delivery of the individual devices remain valid and must be observed by the user of the system.

WARNING: Components and equipment which are not part of this system must not be connected to the system without additional follow-up documentation, assessment and additional notes for the user.

WARNING: The use of endoscopes and light cables which have not been specially approved by KARL STORZ can lead to misdiagnoses.

CAUTION: Only the fluid light cables specially developed for use with the AF, PDD and ICG endoscopes may be used.

Note

The KARL STORZ FI system must not be connected to other devices. All deviations from this instruction require a new assessment of risks and additional documentation with additional warnings and notes for the user if applicable.

All medical electrical components are considered devices for the purposes of EN 60601-1.

All components in the system described are coordinated with one another.

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Safety information for medical electrical systems

A medical electrical system is a combination of several devices of which at least one is a medical electrical device and which are connected via a functional connection, electrical or otherwise, including such connections for transmitting signals and/or substances or via a portable power distributor.

ⓘ NOTE: For the installation of such a system, the safety specifications as per IEC 60601-1-1 or IEC 60601-1 as last amended must be observed.

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Reprocessing

For instructions relating to reprocessing of the system components, please refer to the individual instruction manuals of the system components. Further valid instructions are not intended for the system.

Preventative maintenance measures

Any maintenance work which may be necessary to carry out is detailed in the instruction manual for the individual system components. This includes, for example, the timely replacement of wearing parts. In addition, it is advisable to inspect the supply line of any video cart employed (where applicable) for mechanical damage, and to arrange for a replacement to be carried out by a specialist if damage is detected, before every use.

Once the system has been assembled for the first time with all data lines, signal lines, electrical and non-electrical components, an electrical safety inspection must be carried out. This must be carried out on site or included within the scope of the delivery, if the system has already been inspected for the first time at the factory.

This inspection must be carried out at least once a year. Each alteration to the system (connections, signal lines, cables and components) requires documentation and a new repeat inspection to be carried out according to the following criteria:

These measurements must be carried out as per EN 62353 for type CF devices.

For the permissible limit values for protective ground resistance, earth leakage current, housing leakage current and patient leakage current, please refer to the current standard.

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Maintenance and safety check**Maintenance**

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Safety checks/repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician on this device once a year.

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Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure touch current in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure patient leakage current in accordance with IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

A functional check must be carried out in accordance with the instruction manual.

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during the safety check/repeat inspection, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest edition of the relevant Service Manual.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseas.com

KARL STORZ Endoscopy
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopy Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 Sao Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopy Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4, kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytøften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwon Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411
Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725
Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188
Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renmin Road South
Chengdu, Sichuan, 610016
People's Republic of China
Phone: +86 28 8620-0175
Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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**1 Important information for
users of KARL STORZ
devices**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended as an aid in setting up, connecting and operating the UNIDRIVE® S III motor system with the different application instruments

- Drill Cut-X® II URO handpiece (27702050) with morcellator blade (27056LM)

Please read these instructions carefully. Keep this manual in a convenient and conspicuous location and in its protective casing close to the device.

The instruction manuals/accompanying instruction leaflets enclosed with the respective article must be observed when using the different application instruments.

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

3 Controls, displays,
connectors, and their uses

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- ① Power switch for turning the unit on and off
- ② Bargraph display for indicating the actual speed or maximum speed
- ③ (7-segment) display for
 - maximum speed (specified by operating mode) or
 - error code
- ④ (7-segment) display for
 - selected operating mode
 - set and current speed
- ⑤ Setting the maximum speed of the respective mode
- ⑥ Connection socket for motor/handpiece M
- ⑦ Connection socket for footswitch
- ⑧ Connection socket for footswitch control cord
- ⑨ USB port (for maintenance only)
- ⑩ Potential equalization connection
- ⑪ Power socket
- ⑫ Line fuse holder
- ⑬ SCB connection
- ⑭ Single-pedal footswitch

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Read the instructions carefully before operating
the equipment!

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Connection socket for footswitch/control cord

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Plus and minus buttons for setting the maximum
speed of the respective mode

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Electronic information product pollution control

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DANGER: Risk of explosion if used in the
presence of flammable anesthetic gases.

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CAUTION: Do not open the unit. Danger of
electrical shock!

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Refer to the instruction manual for steriliza-
tion instructions.

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Do not store liquids on or above the unit.

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Keep the unit out of reach of patients.

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8	Configuring the operating mode	22
8.1	DrillCut-X® II URO	22
8.2	Assembly of adaptor.....	24
9	Service and repair	25
9.1	Fuse replacement	25
9.2	Cleaning, disinfection, and sterilization.....	26
9.3	Maintenance and safety check	27
9.3.1	Maintenance.....	27
9.3.2	Safety check.....	27
9.4	Service and repair.....	29
9.5	Disposal.....	29
9.6	Repair program.....	30
9.7	Important information.....	30
9.8	Limitation of liability.....	31
9.9	Warranty policy.....	31
10	Technical description	32
10.1	General.....	32
10.2	UNIDRIVE® troubleshooting table.....	34
10.3	Troubleshooting	36
10.4	Technical data.....	37
10.5	Standard compliance	38
10.6	Directive compliance.....	38
10.7	Technical documentation	39
10.8	Block diagram	40
11	Spare parts, recommended accessories	41
11.1	Spare parts.....	41
11.2	Recommended accessories	42
12	Appendix	43
12.1	Cleaning agents and disinfectants	43
12.1.1	Manual preparation.....	43
12.1.2	Machine preparation	45
12.2	Electromagnetic Compatibility (EMC) information	47
13	Warranty card	63
14	Subsidiaries	66

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5 General information

5.1 System description

The new KARL STORZ motor system UNIDRIVE® S III is characterized by the most modern electronics, flexibility, and efficiency.

The use of separate brushless high-performance electronic motors or brushless high-performance electronic motors integrated in the instrument makes it universally applicable in the application area. The settings and functions appear as 'text' in the display.

The preset parameters are saved, i.e., you only have to select them once. By simply activating the respective menu field, e.g. device information can be displayed. Troubleshooting in the OR is easier since any faults that may occur, such as a motor fault, are displayed automatically.

Special features

- Automatic detection of connected handpieces with program assignment
- Preselectable max. speed
- Constant high motor power over the entire speed range
- Processor-controlled speed and motor torque
- Controls which are clear and easy to read thanks to a color display
- Control output for intended KARL STORZ pump systems

The unit 20701020-1 includes the KARL STORZ-SCB interface (KARL STORZ Communication Bus).

* The SCB interface, based on the CAN field bus, permits remote control of equipment functions and remote display of equipment parameters.

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6.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution, and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

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WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered.

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NOTE: A Note indicates special information about operating the device or clarifies important information.

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WARNING: Read this instruction manual thoroughly before operating the system. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk.

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-  **WARNING:** The UNIDRIVE® S III must only be used in medical facilities that conform to the applicable national regulations.
-  **WARNING:** Keep out of reach of patients.
-  **WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if:
 - they are identified as such in the respective instruction manual; or
 - the intended use and the interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1).
-  **WARNING:** Connect the device to a properly grounded 'Hospital grade' or 'Hospital use' electrical outlet. Check grounds regularly. Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.
-  **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.
-  **WARNING:** Risk of explosion if used in the presence of flammable anesthetic gases.
-  **WARNING:** The UNIDRIVE® S III, in conjunction with the application instruments to be used, is only designed for intermittent operation (50 %, $t=1$ min).
-  **WARNING:** Always unplug the unit from the mains before performing any maintenance or cleaning work on it.

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-  **WARNING:** Do not open device! Risk of electric shock. Refer servicing to the manufacturer or authorized service personnel. Removal of the cover by non-authorized personnel will void the warranty.
-  **WARNING:** Turn the device off if the motor does not start and stop as required.
-  **CAUTION:** The device may only be operated at the voltage stated on the rating plate.
-  **CAUTION:** When replacing fuses, use only fuses of the same rating.
-  **CAUTION:** The UNIDRIVE® S III may only be used with accessories and disposable items which have been designated for use by KARL STORZ.
-  **CAUTION:** Do not store liquids on or above the device. Avoid allowing liquids to enter the unit.
-  **CAUTION:** If the motor becomes blocked for any length of time, it can lead to overheating and cause damage to the motor. If the blockage cannot be disengaged using oscillating operation, the motor must be shut off and the blockage must be removed manually.
-  **CAUTION:** Continuous checks must be performed both before and during use to establish whether the required mode is set.
-  **NOTE:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.
-  **NOTE:** For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices. The current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.

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Before using the UNIDRIVE® S III motor system on a patient, it is imperative that the instructions in this manual are understood and are followed carefully.

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6.2 Intended use

The UNIDRIVE® S III motor system is universally applicable in conjunction with the respective, suitable application instruments in the URO application range:

- During minimally invasive endourological and urological interventions for the controlled excision, morcellation, removal, ablation and transection of tissue in the urinary tract system.
- In urological procedures including nephrectomy.

The exact 'Intended use' of the application instruments is detailed in the respective instruction manual provided with the application instrument.

The UNIDRIVE® S III motor system may only be used with accessories, parts and disposable items that have been designated by KARL STORZ as suitable for the unit.

Unauthorized conversions or modifications to the device are not allowed for safety reasons.

For the intended use of the handpiece 27702050 and the blade 27056LM, please consult the instruction manual enclosed with the handpiece.

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
6.3 User qualification

The UNIDRIVE® S III motor system in conjunction with an application instrument may only be used by physicians and/or medical assistants who have a corresponding specialized qualification and who have been instructed in the use of this device/instrument.

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6.4 Safety precautions at the site of installation

The UNIDRIVE® S III motor system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The UNIDRIVE® S III is equipped with a connector for attaching a ground wire. It should be connected up in accordance with the applicable national regulations.

The UNIDRIVE® S III motor system is not intended for use in hazardous areas. This means, for example:

When using easily combustible and explosive inhalation anesthetics or mixtures thereof, the device must not be operated in the demarcated hazard zone. This also applies for easily combustible and explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

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6.5 Safety precautions when operating the system

WARNING: It is the user's responsibility to make sure the system is safe and operates properly before using it.

During treatment with the UNIDRIVE® S III, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

This device has been tested and is found to comply with the EMC limit values as per EN/IEC 60601-1-2 :2001 CISPR 11 Class B.

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For the safety instructions concerning the handpiece 27702050 and the blade 27056LM, please consult the instruction manual enclosed with the handpiece.

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7 Operating instructions

7.1 Unpacking

Carefully remove the unit and the accessories from the packaging. Check for missing items and for evidence of shipping damage. Please file any complaints with the manufacturer or supplier immediately.

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For use with instrument set

- 1 HOPKINS® 6° wide angle straight forward telescope 27292AMA / 27293AA
- 1 Adaptor 27040LB / 27040SC (for sheath 27050SC)
- 1 Resectoscope sheath 27040SL / 27050SL / 27050SC
- 1 Obturator
- 1 Morcellator blade 27056LM

① NOTE: Optional accessories can be found under 'Recommended Accessories'.

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7.3 Installing and connecting the unit


NOTE: The UNIDRIVE® S III including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

WARNING: The UNIDRIVE® S III is not intended for use in potentially explosive atmospheres. Do not operate the system within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Place the UNIDRIVE® S III on a level surface. The UNIDRIVE® S III is equipped with a connector for attaching a ground wire (14). The unit's ground line should be installed by a qualified electrician.

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Before connecting the UNIDRIVE® S III to the electrical supply line, verify that the line voltage stated on the unit's identification plate is the same as that of the electrical supply line to be used.

NOTE: At the factory, 2 A fuses suitable for operation at a 230/240 V supply are inserted in the unit.

If you are operating the device with a 100...120 V power supply, the 2 A fuses inserted in the factory must be exchanged for 4 A fuses (T4AH or T4AL) (see Section 'Fuse replacement', p. 25).

Connect the power, insert the power cord end into the receptacle (15) of the unit.

WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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
7.4 Connecting the footswitch


To use the footswitch ⑪, insert the connecting cord* of the footswitch into the socket for the footswitch ⑫ on the front of the UNIDRIVE® S III. Align the red dot on the connecting cord with the red dot on the footswitch socket.

NOTE: If no footswitch is connected when the unit is switched on, the message: 'E11' flashes in the 7-segment display ④.

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7.5 Connecting application instruments

Prepare the application instrument according to the instruction manual.

NOTE: Read the accompanying instruction leaflet provided with the application instrument carefully!


Connect the connecting cable* coming from the application instrument (morcellator handpiece) to the motor output ⑥.

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The connection will be made properly when the red dot on the connection cable is aligned with that on the socket.


NOTE: If no motor is connected when the unit is switched on, the message: 'E12' flashes in the 7-segment display ④.

*** NOTE:** To unplug the connecting cable, pull out the plug — do not pull directly on the cable!

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
7.6 Combined operation with ENDOMAT® LC

- ① NOTE:** The respective instruction manual for the combi unit must be read and followed extremely carefully prior to the use and connection of the combi unit on the motor control unit UNIDRIVE® S III and connected accessories.
- ① NOTE:** The combi unit ENDOMAT® LC must only be used as a suction pump with the suction tube set 20330343.

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Connect the ENDOMAT® LC to the combi unit using the footswitch control cord 20701070.

To do so, insert the footswitch control cord into the pump connector (⑧) socket on the rear of the UNIDRIVE® S III.

Insert the other end of the control cord on the back of the unit into the footswitch input socket (see figure) of the ENDOMAT® LC.

Align the red dot on the connecting cord with the red dot on the footswitch socket and push in as far as it will go.

The ENDOMAT® LC must only be operated via the footswitch control cord and the footswitch connected to the UNIDRIVE® S III.

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- ① Pump tube
- ② Narrow tube connector
- ③ Wide tube connector
- ④ Suction line (part 2) to suction bottle
- ⑤ Suction bottle connection
- ⑥ Suction line (part 1) to handpiece

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Connect both ends of the pump tube to the connectors until they click into place. Connect the inflow tube coming from the DrillCut-X® II to the narrow pump tube connector. Attach outflow tubing to the wide pump tubing connector.

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- ① **NOTE:** If you are using a collecting bottle without KARL STORZ connections or a tissue collection tray, the connectors ⑤ must be removed.
- ① **NOTE:** Observe 'Leakage testing' on page 37 of the ENDOMAT® LC instruction manual.

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- ① **NOTE:** The tubing sets 20330343 are delivered unsterilized. For sterilization, see page 26 ff.
- ① **NOTE:** The tubing is delivered ready-assembled. If the tubing has been disassembled, for preparation for example, it must be reassembled as shown in the figure.

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7.6.4 Inserting the pump tube

Flip the lever for retracting/extending the pump rollers ⑧ down into the horizontal position.

CAUTION: Only use tubing sets designed for the ENDOMAT® LC.

CAUTION: Make sure that the tube is completely inserted into the pump bed.

NOTE: The pump tubing connectors have different widths so that the pump tubing can only be inserted in one direction.

Attach the pump tubing to the roller pump and press both connectors into the fastening grooves.

WARNING: Due to application wear, we recommend changing the pump tubing after 5 uses.
(Order no. 20.3303.93, pack of 25)

CAUTION: If the pump tubing is perforated, there is a danger of fluid entering the device!

Press the lever for retracting/extending the pump rollers ⑧ up into the vertical position; the pump tubing is locked in place.

NOTE: If the pump tube is not locked in place, the pump cannot be started.

7.6.5 Removing the pump tube

Press lever ⑧ down into the horizontal position.
Remove the pump tube from the reel pump.

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7.6.6 Initial operation

Pump lever closed.
Switch on the power switch ①.

WARNING: The user must monitor the display test when turning on, to ensure that all displaying elements are in perfect working order.

NOTE: Once the self test at power up has successfully completed, an acoustic availability signal sounds.

The setpoint for the flow flashes in the display (digital display ②).

NOTE: Once switched on, the unit awaits confirmation or a change of the last setpoint used. This is shown flashing. The LED in the irrigation/suction START/STOP button ⑥ also flashes. The pump cannot be started as long as the displays are still flashing.

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7.7 Setpoint confirmation


To confirm the setpoint, press the START/STOP button ⑥.

The LED in the START/STOP button ⑥ goes out and the digital display ② shows the set point continuously.

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7.8 ENDOMAT® LC operating modes

The ENDOMAT® LC distinguishes between two different modes:

- normal mode and
- high flow mode (bolus function)

You can change between the modes to enter the setpoints via the M button ⑤.

NOTE: The high flow function can only be activated via the corresponding external switch on socket ⑨.

NOTE: During combined operation with the UNIDRIVE® S III, high flow mode is not possible.

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7.8.1 Select normal mode flow rate

In normal mode (without symbol display) it is possible to select a flow rate of 0-1,000 ml/min (in increments of 50 ml).

Press the ± buttons ③ until the digital display ② displays the required setpoint.

NOTE: If the ± buttons are held pressed for longer than 3 s, the values continue to change until released.

NOTE: A setting of 1,000 ml/min is recommended when using morcellator blades.

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
The pump starts in normal mode

The pump can be started using the single-pedal footswitch 20016230 (1st step) connected to the UNIDRIVE® S III.

NOTE: Section 7.9 onwards of the instruction manual for the UNIDRIVE® S III must be read prior to activating the footswitch.

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7.8.2 Test for proper operation and air-relief of the tubing system

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.

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7.8.3 Checking the suction

Close off the suction tube (instrument tube) with your finger and check whether a vacuum forms.

WARNING: If the reaction of the device deviates significantly, it must be checked by an authorized technician before being used again.

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7.9 Initial operation


Switch the UNIDRIVE® S III on at the power supply ① (power switch is illuminated).

- ① **NOTE:** Before switching the unit on, the user must connect the footswitch and an application instrument to the unit. The software is designed to verify the connections and will indicate an error message in the 7-segment display ④ if necessary.

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
The screen shown on the left appears for the time it takes the application to start. A self-test is then performed in which all the LEDs illuminate for a short time.

- ① **NOTE:** The next pages contain a description of the simplest application example; adopting the most recently used mode, as well as the adjustments for speed. During the switch-on and switch-off procedures, a single signal tone is heard.

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7.10 UNIDRIVE® S III function

A start screen is displayed once the unit has started. The function only becomes active once application instruments have been connected with the unit.

The operation screen now shows the following factory settings: Operating mode 'Morcellator', motor active, oscillation. This operating mode is displayed in the 7-segment display ④ by P6.

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
7.11 Motor connection

Various application instruments operated within various speed ranges can be connected up to the control unit. To do so, the required instrument must be connected to the motor connection. Having detected the handpiece automatically, the UNIDRIVE® S III loads the standard parameters and displays them in the 7-segment display ⑤.

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7.12 Overview of application instruments

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Application instrument	Program	Operating mode	Speed range
DrillCut-X® II URO handpiece (27702050) with morcellator blade (27056LM)	P6	DrillCut-X® II URO	500 – 5,000 rpm Preset = 4,000 rpm

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8 Configuring the operating mode

8.1 DrillCut-X® II URO

After switching on the UNIDRIVE® S III, the adjacent screen appears. The handpiece connection cable already inserted in the device is detected and the operating mode Morcellator P6 automatically assigned and activated.

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Configuring the operating mode: Preselecting the maximum speed (setpoint value)

Use the '+/- buttons' ⑥ to define the maximum speed (set value) (display ⑤). The available speed range is displayed in the bargraph display ②. The preset range is displayed by an LED in the bargraph display ②.

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If the 'handpiece' is activated by pressing the foot-switch pedal to the second position, the bargraph display ② and the display ⑤ are updated to the preset value.

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Footswitch function

The footswitch has the following functions in the 'URO' operating mode:

The footpedal has a two-step function. When the pedal is pressed up to the point of resistance (approx. up to 3/4 of the pedal movement range), only the suction pump connected externally via a connecting cable is activated.

When the point of resistance is passed, the morcellator function is also activated in addition to the suction pump in the last 1/4 of the pedal movement range.

The separate suction function of the first position allows targeted suction of the tissue to be extracted prior to activation of the cutting mode and targeted suction of loose tissue concretions without cutting mode.

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
8.2 Assembly of adaptor

1. Assembly of adaptor 27040LB/27040SC on telescope 27292AMA/27293AA
 2. Assembly of telescope with adaptor on sheath (27040SL/27050SL/27050SC)
 3. Assembly of irrigation tube on telescope/sheath
 4. Assembly of suction tube 20330343 on DrillCut-X® II URO handpiece 27702050
 5. Assembly of morcellator blade 27056LM on handpiece
 6. Insertion of blade in the telescope's working channel
 7. Connection of suction tube to a collection system
- Connection of the optional components in suction mode
to the mtp* tissue collection tray +
mtp* container system, as well as
KST reusable 5 l suction bottle
20300050 + cap 20300034 +
mtp* filter

NOTE: Ensure the possibility of venting the container system.

NOTE: The instruction manual enclosed with the morcellator handpiece 27702050 must be followed extremely closely!

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Remove the power fuse holder ⑮ with a screwdriver or another suitable tool.

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! CAUTION: Use only fuses with the correct nominal ratings.


Insert new fuses.
Insert the power fuse holder again.

	100...240 V
Line fuse	2xT 4.0 AH 250 V

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Reconnect the power.
Check that the instrument works properly.

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9.2 Cleaning, disinfection, and sterilization

-  **NOTE:** Please read the instruction manual provided with the application instrument carefully!
-  **NOTE:** The 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions must be followed. The cleaning, disinfection, and sterilization procedures are explained therein in detail.
-  **NOTE:** The instructions for cleaning and sterilizing the combi device used in combination with the UNIDRIVE® S III and its accessories must be taken from the respective instruction manual for the combi device.
-  **NOTE:** For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices. The current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.
-  **WARNING:** Always unplug power before cleaning.
-  **CAUTION:** Avoid allowing liquids to enter the housing at all costs.

Wipe-clean the exterior surfaces of the unit with a disposable cloth moistened with disinfectant. For wipe-down disinfection, we recommend using only surface disinfectants (no solutions predominantly containing alcohol such as rapid-action disinfectants) and 'disposable disinfectant methods' based on ammonium chloride.

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9.3 Maintenance and safety check

9.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

9.3.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

8. Check the device and accessories for any mechanical damage which may impair functionality.
9. Check that inscriptions relevant to safety are legible.

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Electric measurements:

- Inspect the device safety fuses
- Protective ground resistance measured according to IEC 60601-1/A2:1995: Limit value: 0.2 ohms
- Earth leakage current measured according to IEC 60601-1/A2:1995: Limit value: NC: 500 µA, SFC: 1000 µA
- Patient leakage current measured according to IEC 60601-1/A2:1995: Limit value: 100 µA (BF), 10 µA (CF)
- Leakage currents must not exceed the limit values specified above.

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The safety check must be entered in the unit log-book and the results of the check must be documented.

If the device does not function reliably and/or safely, it must be repaired.

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9.4 Service and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine KARL STORZ parts only.

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9.5 Disposal

This device has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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9.6 Repair program

During the repair period, you will generally be loaned a unit, which must then be returned to KARL STORZ as soon as the repaired unit is returned to you.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the unit has been used in accordance with its operating instructions at all times.

9.9 Warranty policy

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot bear the cost of transportation or the risks associated with shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the warranty card and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

NOTE: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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10 Technical description

10.1 General

In the event of the unit malfunctioning, the trouble-shooting list provided in this chapter will assist the user in finding the possible cause. The trouble-shooting list indicates how to find and rectify the various malfunctions.

WARNING: High voltage is present in the unit's interior. The unit must only be opened by an authorized KARL STORZ Customer Service technician.

WARNING: Maintenance work that is not mentioned in this chapter must only be performed by technical personnel that are authorized by KARL STORZ.

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General information

The UNIDRIVE® S III displays error messages in the 7-segment display ④. Hand the device in for servicing when a hardware error message is displayed (E01-E10). Take appropriate action for user error messages (E11-E30). E.g. E11 'NO FOOTSWITCH' – check for whether the footswitch is properly connected. Detailed information can be obtained from the error table on page 36 (also see the 'Service mode' section).

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10.2 UNIDRIVE® troubleshooting table

① NOTE: If the device is not working as required, check the error message displayed in the 7-segment display. If necessary, contact the KARL STORZ Service.

Error message	Cause	Remedy
1 = E01: MEMORY ERROR	Parameters could not be saved successfully.	Turn device on/off. Load default values. Or send unit in for repair.
2 = E02: MOTORBOARD DEFECTIVE	Motor control overloaded.	Turn device on/off. Or send unit in for repair.
3 = E03: HARDWARE CONFIGURATION ERROR	Hardware configuration is incorrect	Send unit in for repair.
4 = E04: SOFTWARE CONFIGURATION ERROR	Software configuration is incorrect.	Send unit in for repair.
5 = E05: –		
6 = E06: –		
7 = E07: –		
8 = E08: –		
9 = E09: MOTOR DEFECTIVE M	Motor on channel 1 is defective.	Release footswitch and restart. If the error persists, a different handpiece must be used.
10 = E10: –		
11 = E11: NO FOOTSWITCH	Footswitch not connected correctly.	Connect footswitch. If the error persists, check whether the error is on the footswitch or unit. Send the defective component in for repair.
12 = E12: NO MOTOR	No motor/handpiece connected.	Connect motor/handpiece.
13 = E13: TOO HIGH TEMPERATURE MOTOR M	Motor/handpiece on channel 1 overheated.	Allow the connected handpiece to cool down or operate an additional handpiece on channel 2.
14 = E14: –		
15 = E15: –		
16 = E16: –		
17 = E17: MOTOR NOT SUPPORTED M	The motor/handpiece connected on channel 1 is not supported.	Only use handpieces approved by KARL STORZ. A service update may be necessary.
18 = E18: –		
19 = E19: MOTOR STOPPED M	Motor/handpiece on channel 1 blocked in operation.	Release footswitch and restart. If the error persists, clean the handpiece or exchange the attachment.
20 = E20: –		

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Power switch lamp does not light up.

Possible cause

- 1 Power supply failure.
- 2 Defective power fuse.
- 3 Power switch lamp defective.

Recommended remedy

- 1 Check that there is electricity to the wall outlet.
- 2 Change fuses as described in the instruction manual. Make sure you fit the correct type of fuse.
- 3 Notify customer service.

Problem

- 1 Power switch lamp is illuminated.
- 2 No indication of rotation.
- 3 No energy output when operating the foot-switch/buttons.
- 4 Speed cannot be regulated.

Possible cause

- 1 Faulty motor control electronics.
- 2 Connection to motor interrupted.
- 3a Faulty handpiece.
- 3b Handpiece blocked or jammed.
- 4 Faulty motor control electronics.

Recommended remedy

- 1, 2, 4 Notify customer service.
- 3a Connect another handpiece.
- 3b Clean the handpiece.

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10.4 Technical data

UNIDRIVE® S III	20701020-1
Line voltage	100...120/ 230...240 V~
Line frequency	50/60 Hz
Line fuse	2xT 4.0 AH 250V, IEC 127 format
Power consumption	max. 200 VA
Operating temperature	10-40 °C
Humidity	30 % - 75 %
Dimensions	305 x 165 x 233 mm (W x H x D)
Weight	4 kg
Operating mode (in connection with application instrument)	intermittent periodic duty
Duty cycle	50 %; t = 1 min
Storage/transport conditions:	
Humidity (non-condensing)	20 %...95 %
Temperature	(-)10 °C...60 °C
Atmosph. Pressure	500hPa...1080 hPa
DrillCut-X II URO	
Minimum speed:	500 rpm
Maximum speed:	5,000 rpm
Rotation mode:	Oscillation

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(for 20701020-1)

According to IEC 60601-1, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shock: Protection Class I
- Degree of protection against electric shock: Applied part type BF
- Type of protection against moisture: Drip-water protection as per IPX 1

Conforms with IEC 60601-1-2 :2001:

Observe the information on electromagnetic compatibility in the appendix (pages 47-59).

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10.6 Directive compliance

Medical product of class IIa

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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**11 Spare parts,
recommended accessories**

11.1 Spare parts

Part	Order no.
Tube connectors to suction tube (203300343) to suction bottle	20300180
Tube connector set to suction tube (203300343) to reel pump	20300482
Replacement pump tubes to suction tube (203300343), 25 pieces	20330393
Replacement ball for overflow protection device on bottle lid 20300034, two required	2911590
Replacement sleeve for overflow protection device on bottle lid 20300034,	20300580
Line fuse T 4.0 OAL 250 V (100...240V~)	2027690
Mains cord (grounded)	400 A
'Hospital Grade' mains cord (USA)	400 B
Instruction manual	96156019D

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11.2 Recommended accessories

(Not included in the scope of supply)

Part	Order no.
Universal spray, pack of 6 cans	280052
Suction tube to ENDOMAT LC	20330343
Single-pedal footswitch, two-step, with switching function	20016230
SCB connecting cable	20090170
Footswitch control cord	20701070
HOPKINS 6° wide angle straight forward telescope	27293AA
Resectoscope sheath	27050SL
	27050SC
Adaptor	27040LB
	27040SC
Obturator	27040OC
Morcellator blade	27056LM
mtp suction container	mtp*
5 l suction container	20300050
Suction bottle lid	20300034
DrillCut-X® II handpiece	27702050
Disposable tissue collection tray, pack of 10	030970-10*

*mtp medical technical promotion gmbh,
Take-Off GewerbePark 46
78579 Neuhausen ob Eck, Germany
Phone: +49 (0)7467 945 04-0
Fax: +49 (0)7467 945 04-99

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12 Appendix**12.1 Cleaning agents and disinfectants**

Chemicals for the preparation of telescopes and endoscopic instrument sets

12.1.1 Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid:

- HAMOU telescopes with eyepiece drive mechanism
- Magnifiers with drive mechanism
- IMPERATOR drill handpieces
- Prismatic light deflectors

① NOTE: For motors and their handpieces, the instructions in their respective instruction manuals must be followed.

The following instruments produced by KARL STORZ are not suitable for cleaning in an ultrasound cleaning bath:

- Rigid telescopes
- Flexible telescopes
- Other instruments with telescopic glass components

When preparing and using the solutions, follow the chemical manufacturers' instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer.

① NOTE: For videoscopes please refer to our separate list of approved chemicals.

The materials approved by KARL STORZ as regards material compatibility can be found on the following page.

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Typ

- 1 = Cleaner
2 = Disinfectant
3 = Suitable for ultrasound baths
A = Not suitable for flexible endoscopes
(A) = Only for new generation flexible endoscopes (year of manufacture >2005)

NOTE: You will also find the current list of approved materials on the Internet at www.karlstorz.com.

NOTE: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

CAUTION: The use of CAUSTIC SODA can cause changes to the surface of aluminum parts (including coated parts), plastics, and soldered connections, and may reduce the service life of the instrument set.

For flexible endoscopes, only immerse the examination shaft, not the housing/handle part.

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12.1.2 Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a washer and disinfecter:

- HAMOU telescopes with eyepiece drive mechanism
- Magnifiers with drive mechanism
- IMPERATOR drill handpieces
- Standard, tungsten carbide, and diamond drills
- Prismatic light deflectors

ⓘ NOTE: For motors and their handpieces, the instructions in their respective instruction manuals must be followed.

Selection of a method for cleaning and disinfection of telescopes must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special procedures must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable.

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12.2 Electromagnetic Compatibility (EMC) information

WARNING: Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The UNIDRIVE® S III model 20701020-1 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The UNIDRIVE® S III model 20701020-1 is a Group 1 unit (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device;
- increase the distance between the equipment;
- connect the devices to different electrical circuits.

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WARNING: The UNIDRIVE® S III model 20701020-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

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WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the UNIDRIVE® S III model 20701020-1. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 201		
Guidance and manufacturer's declaration – electromagnetic emissions		
UNIDRIVE® S III model 20 701020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III model 20 701020-1 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	UNIDRIVE® S III model 20 701020-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The UNIDRIVE® S III model 20 701020-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class D (see note below)	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	not applied (see note below)	
 NOTE: The measurement of harmonic emissions is not required (power < 75 W). The measurement of voltage fluctuations/flicker is not required (activated manually, unlikely to produce significant voltage fluctuations or flicker).		

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Table 202			
Guidance and manufacturer's declaration – electromagnetic immunity			
UNIDRIVE® S III model 20701020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III model 20701020-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycle 70 % U_T (30 % dip in U_T) for 25 cycle <5 % U_T (>95 % dip in U_T) for 5 s	0 V, 10 msec 92 V, 100 msec 168 V, 500 msec 0 V, 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Table 204 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
UNIDRIVE® S III model 20701020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III model 20701020-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF disturbances acc. to IEC 61000-4-6 Conducted RF disturbances acc. to IEC 61000-4-3	3 V_{ms} 150 kHz to < 80 MHz 3 V/m 800 MHz to 2,5 GHz	3 V_{ms} 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the UNIDRIVE® S III model 20701020-1/ 40701720-1/40701820-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz $d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [3.5/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 206
Recommended separation distances between portable and mobile RF communications equipment
and the UNIDRIVE® S III model 20701020-1

UNIDRIVE® S III model 20701020-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \sqrt{\frac{P}{f}}$	80 MHz to < 800 MHz $d = \sqrt{\frac{P}{f}}$	800 MHz to 2.5 GHz $d = \sqrt{\frac{P}{f}}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

① NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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x

WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

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To be completed by the supplier/importer:
Company stamp/signature:

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karl-storz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseaa.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karl-storzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseaa.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América do Sul Ltda.
Rua Joaquim Floriano, 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karl-storz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karl-storz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsensvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karl-storz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skarholmen, Sweden
Postal address: P.O. Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karl-storz.se

KARL STORZ Endoscopy Suomi Oy
Valimotie 1 B. 4, kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karl-storz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlystvej 33
2840 Høfte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karl-storz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karl-storz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopier.be
KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karl-storz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karl-storz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karl-storz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karl-storz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karl-storz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia

Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karl-storz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colanin, nr 74, Sector 4
041393 Bucuresti, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karl-storz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Gounan 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karl-storz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy
000 KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karl-storz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ Endoskope
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391
E-Mail: info@karl-storz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karl-storz.co.za

TOO KARL STORZ Endoscopy Kazakhstan
Sary-Arka str., 6, BC "Arman", off.430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karl-storz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karl-storz-emg.com

KARL STORZ Endoscopy India Private Ltd
D-181, Oldia Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karl-storz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F10-Q, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karl-storz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1, Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karl-storz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
2012 3301 Beijing, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renmin Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd,
Shenyang Branch Company.
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karl-storz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karl-storz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jelton Tower 1, Suite 905
16 Neudi-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.lee@karl-storz.com

KARL STORZ Endoscopy Taiwan Ltd
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karl-storz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karl-storz.co.jp

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
Important information for users of KARL STORZ devices

Thank you for your confidence in the brand name KARL STORZ. Like all of our products, this product too is the result of years of experience and careful manufacturing. You and your organization have decided in favor of a modern high quality product from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the EndoCrane™ holding arm system. All the necessary details and actions are clearly explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

1. General warning

- The product must only be operated as intended and in accordance with the instruction manual and only together with correspondingly trained and qualified specialist staff!
- Maintenance and repairs must only be performed by authorized specialists. For more information, please contact the respective representative.
- The product can only be operated in combination with the accessory components specified in the instruction manual. Other combinations and accessory components may only be used if they are expressly suited to the intended application and do not affect performance characteristics or safety requirements.
- The products must be treated as described in the instruction manual prior to every use and shipping for the protection of the patient, user, and third parties.

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• To avoid damage: if the holding arm system is not required, it should be decommissioned and kept movable (photos see p. 24)

- Switch off ON/OFF rocker switch (17)
- Remove connector (15) from the connection socket (14)
- Connect connector (15) to emergency release plug (18)
- Press release button (6)

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

- 1 Base 28172 HG, with 2 thumb screws, for fastening on operating table
- 2 Stand bar
- 3 Union nut
- 4 Holding arm (already assembled when delivered, cannot be disassembled)
④ + ② = 28272ECS
- 5 Handle
- 6 Release button (must be activated for the holding arm to be movable) *
- 7 Joint 1 (not visible in photo)
- 8 Joint 2
- 9 Joint 3
- 10 KSLock (patented quick-release coupling)
- 11 Clamping jaw (here 28272UGN; alternatives: clamping jaws 28272 UG / UK / UU). The clamping jaws are not supplied as standard.
- 12 Standard device rails
- 13 Control unit (also called holding arm control unit)
- 14 Connection socket
- 15 Connector *
- 16 Supply cable
- 17 ON/OFF rocker switch *
- 18 Emergency release plug *
- 19 Supply cable connection; cannot be disassembled
- 20 Fastening screw for base
- 21 Stand bend

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- ⑰* ON/OFF rocker switch
- ⑳* Alarm silent (not illuminated)
- ㉑ Unit ON
- ㉒ Observe the instruction manual
- ㉓ Standby mode ON / OFF (not illuminated)
- ㉔ Standby mode ON / OFF (illuminated)
- ㉕ Alarm display (illuminated)
- ㉖ Holding arm symbol (with illuminated dot)
- ㉗ Connection socket
- ㉘ Unit type BF
- ㉙ Potential equalization
- ㉚ Power switch
- ㉛ Fuse holder
- ㉜ Power supply

* For positions ⑱ to ㉑, please see previous page.

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Read the instructions carefully before operating the
equipment

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Potential equalization connector

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1. General warning	IV
2. Photographs of the unit	V
3. Safety instructions	3
3.1 Warnings and cautions	3
3.2 Definitions	3
3.3 Warnings	4
3.4 Intended use	9
3.5 User qualification	9
3.6 Possible combinations	9
3.7 Safety precautions at the site of installation	10
3.8 Safety precautions when operating the device	10
4. Installation and operating instructions	11
4.1 Unpacking	11
4.2 Basic equipment	11
4.3 Assembly of the EndoCrane™ holding arm system and commissioning	12
4.3.1 Visual inspection of the individual components	14
4.3.2 Perform the following visual inspections:	14
4.3.3 Construction of the holding arm system	14
4.3.4 Construction of the stand	15
4.3.5 Assembling the holding arm	16
4.3.6 Connecting the holding arm to the holding arm control unit	16
4.3.7 Checking the holding force	19
4.3.8 Functions of the handle	20
4.3.9 Sterile covering of the holding arm	21
4.3.10 Fastening the clamping jaw	22
4.4 Disassembly of the EndoCrane™ holding arm system	24
4.5 Important notes on avoiding damage	24
4.6 In case of fault	25
4.6.1 Consequence of power failure	26
4.7 Procedure in case of fault	27
5. Alarm messages, troubleshooting and remedies	28
5.1 Alarm messages	28
5.2 Troubleshooting and remedies	28
5.3 Storage when normal operation is interrupted	29
5.4 Replacing spare parts	29

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6. Cleaning, disinfection, sterilization	30
6.1 Preparation of the holding arm control unit by the user	30
6.2 Preparation of the holding arm system by the user	30
6.2.1 Removal of the sterile drape	30
6.2.2 Disinfection of the holding arm and stand	31
6.2.3 Steam sterilization (not for the holding arm!)	31
6.2.4 Preparation of the clamping jaws and the base	32
7. Service and repair	33
7.1 Maintenance and safety check	33
7.1.1 Maintenance	33
7.1.2 Safety check	33
7.2 Disposal	35
7.3 Repair program	35
7.4 Important information	35
7.5 Limitation of liability	36
7.6 Warranty policy	36
8. Technical description	37
8.1 Technical data for the holding arm control unit	37
8.2 Interfaces to the holding arm control unit	37
8.3 Technical data for the holding arm with stand	38
8.4 Operating, storage, and transport conditions	38
8.5 Standard compliance (for 20 7800 20)	39
8.6 Directive compliance (for 20 7800 20)	39
8.7 Technical documentation	39
8.8 Software user rights	40
9. Appendix	41
9.1 Electromagnetic compatibility (EMC) information	41
10. Subsidiaries	56

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3.1 Warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words Warning, Caution, and Note convey special meanings. Wherever they are used in the instruction manual, please read the text that follows carefully, in order to ensure safe and efficient operation of the device.

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3.2 Definitions

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: Notes contain additional information on the safe operation of the device.

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WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using the EndoCrane™ holding system. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk. Prior to the first application of the EndoCrane™ holding system on the patient, it is imperative that you acquaint yourself with the functioning and operation of the system. The instruction manual should be made available to every authorized person of the operating personnel who is qualified to use the EndoCrane™ holding arm system.

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3.3 Warnings

WARNING: The instruction manuals and interface specifications for the medical devices used in combination must be observed precisely. Electromagnetic and other interferences may arise between the product and other products, which could cause disruptions and/or malfunctions.

WARNING: Test the holding arm system before each use to ensure that it functions correctly.

WARNING: Risk of explosion if the unit (control unit ) is used in the immediate vicinity of a flammable atmosphere (e.g., anesthetics). The unit is not explosion-protected.

WARNING: Danger of an electric shock! Do not open the control unit ! Refer servicing to authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims.

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WARNING: Grounding reliability can only be achieved when the control unit  is connected to a Hospital Only  Hospital Grade outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: The electrical installations in the operating room in which the unit control is connected and operated must comply with the applicable IEC standards.

WARNING: Prior to all work on the control unit  (e.g., cleaning, servicing), the unit must be switched off and disconnected from the power supply.

WARNING: Keep control unit  out of reach of patients.

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Safety instructions Warnings and cautions

WARNING: The device may only be operated at the voltage stated on the rating plate.

WARNING: When replacing fuses, use only fuses of the correct rating. The fuse must only be replaced by qualified personnel.

WARNING: The unit must only be connected to the supplied power cord or a power cord with the same specifications.

WARNING: Prior to every use, the holding system should be checked for defective rough surfaces, sharp edges, and projecting parts in order to rule out safety risks.

WARNING: Prior to every use, the units, instruments, and accessories should be checked for damage, soiling, and completeness. All connecting cords must be checked for damage. Any inscriptions must be complete and easily legible.

WARNING: Risk of infection for patient/personnel: These instruments are **not sterile** when delivered. Clean, disinfect, and sterilize* the instruments before initial use and before every subsequent use using validated procedures.

* Only the base ① and the clamping jaw ⑪ may be sterilized!!

WARNING: The instruction manual must always be observed when cleaning the holding system and preparing/sterilizing the unit accessories.

WARNING: As the sterile drape ⑫ for the holding arm is usually a disposable item, it must be disposed of in accordance with the applicable regulations for the disposal of medical disposables. Further information can be found in the current specialist literature. Observe the drape manufacturer's specifications.

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WARNING: Photo 45: Changing the position of the base ① and the stand ② also changes the position of the instrument/endoscope and/or endoscopic accessory. This could result in injury to the patient! When positioning the holding arm above the patient, the fastening screws ② on the base and the union nut ③ on the stand must not be loosened!

WARNING: Photo 46: Test that the lock is functioning correctly before every use: Once the supply electronics are turned on, the holding arm should lock. Press the two release buttons ⑥ located on opposite sides of the handle ⑤ several times and move the holding arm ④. Repeat the procedure several times. While doing so, check that the holding arm is functioning perfectly (release/lock).

WARNING: Photo 47: When releasing the arm via the two release buttons ⑥, it must be ensured that the holding arm ④ does not move in an undesired or uncontrolled manner! The user must support the holding arm before pressing the release buttons, if this is necessary.
After the release buttons have been released, it must be checked that the arm locks properly!

WARNING: In case of system failures: the emergency lock will last for at least 10 minutes! After that, the holding arm may drop! (See arrow)
Decommission the EndoCrane™ holding arm system immediately if system failures occur and remove it from the operating field (see page 29).

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
WARNING: Photo 48: Releasing the clamping jaw can alter the position of the instrument/endoscope and/or endoscopic accessory. This could result in injury to the patient!

WARNING: Never perform the following procedure while the holding arm system is above the patient! **Photo 49:** Before you loosen the fastening screws of the base in order to reposition the holding arm system, you must support the holding arm's center of gravity so that once the fastening screws have been loosened, you can take the whole weight of the holding arm in your hand.

WARNING: Photo 50: The joint locks can be completely released by connecting the emergency release plug to the connector and then pressing the two release buttons on the handle at the same time!

When the holding arm is uncoupled to avoid damage and for better handling, always connect the emergency release plug to the connector!

WARNING: Photo 51: All the elements of the holding system which can be tightened (e.g., levers, screws), should only be tightened manually (do not use any tools). Otherwise, the elements of the system may be damaged, which could in turn endanger the patient or it make it impossible to loosen the elements.

WARNING: Do not bend the cable too hard. Store the cable in loops at least 15 cm in size. Cables stored on the floor must not be stepped on, squashed, or rolled over with objects. Risk of cable breaking!

WARNING: Please observe the instruction manual for the base fastened on the operating table.

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Safety instructions
Warnings and cautions

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CAUTION: The system's connection cord should always be layed in such a way that it cannot be damaged by the personnel or instrument set.

CAUTION: Do not prepare the holding arm in steam and/or hot air sterilizers. Prior to use in the OR, cover with a sterile drape. Only the clamping jaw and the base supplied for fastening to the operating table may be treated with the steam sterilization procedure.

CAUTION: No cleaning agents, abrasives, or solvents may be used in the care of the unit (control unit). Observe the manufacturer's specifications.

CAUTION: To prevent damage to the products during transportation, we recommend using the original packaging for shipping. Always equip the holding arm with the emergency release plug and release it via the button on the handle before shipping.

NOTE: Please specify the type number and serial number provided on the manufacturer's identification plate with all inquiries and in all written correspondence.

NOTE: Every time the holding system is commissioned, the steps must be performed as carefully as they are for the initial commissioning. The listed assembly steps, warnings, cautions, and notes must be observed for every commissioning.

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Before using the unit on the patient, it is imperative that you acquaint yourself with how the unit operates and is controlled.

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3.4 Intended use

The EndoCrane™ holding system was developed for endoscopic interventions and is used for the holding and one-handed positioning of endoscopes and endoscopic cameras in diagnostic and treatment applications.

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Contraindications: No contraindications relating to the product are currently known. Use of the instrument is contraindicated if, in the opinion of the responsible physician, such use would endanger the patient. E.g., due to the general condition of the patient, or if the endoscopic method itself is contraindicated.

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3.5 User qualification

The EndoCrane™ holding arm system may only be used by physicians and medical assistants with the appropriate specialized qualifications and who have been instructed in use of the unit.

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
3.7 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

It is not intended for use in hazardous zones. This means, among other things:

Wherever easily combustible and explosive inhalation anesthetics and their mixtures are used, the device may not be operated in the demarcated hazard zone. This also applies to easily combustible and explosive chemicals, e.g. skin disinfectants and surface disinfectants.

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3.8 Safety precautions when operating the device

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

During treatment using the EndoCrane™ holding system, the patient must be treated and kept under observation with the usual medical care. This includes monitoring sterile application conditions where required by the type of intervention.

This device has been tested and complies with the EMC limits values as per EN 60601-1-2 / IEC 601-1-2.

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Tel.: ++49 7467 945040
Fax: ++49 7467 9450499
URL: www.mtp-tut.de

**4. Installation and operating
instructions**

4.1 Unpacking

Carefully unpack the EndoCrane™ holding arm system and remove the unit and its accessories from their packaging. Check the delivery for missing items and possible damage.

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4.2 Basic equipment

Control unit	20780020
Power cord	400A (240V) Intl. 400B (110V~) U.S.
Stand and holding arm	28272 ECS
Base, for fastening on operating table	28172 HG
Spring balance	28272 ECW
Drape, sterile, pack of 25 pieces	040146-25*
Instruction manual	96116038D

* mtp medical technical promotion gmbh
Take-Off GewerbePark 46,
78579 Neuhausen ob Eck
Germany

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**4.3 Assembly of the EndoCrane™
holding arm system and
commissioning**


NOTE: The EndoCrane™ holding arm system may only be used in medical rooms where all electrical equipment has been installed in accordance with applicable national regulations.

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WARNING: The EndoCrane™ holding arm system is not intended for use in potentially explosive atmospheres. Do not operate the unit within the demarcated hazard zone while explosive anesthetic gases are in use.

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Initial assembly and commissioning

There is no difference between the initial assembly and commissioning and subsequent commissioning. This means that the following points must be observed every time.

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WARNING: Every commissioning must be performed in the same way. For this reason: the listed assembly steps and warnings, precautions, and notes must be observed for every commissioning.

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
Place the holding arm control unit on a flat surface.

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
The holding arm control unit is equipped with a connector for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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The device may only be operated at the voltage stated on the rating plate.

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Connect the power cord and insert the power plug fully into the power socket ④.

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4.3.1 Visual inspection of the individual components

Prior to every assembly of the holding arm, all components should be inspected visually so that damage or soiling is identified in advance. A test for proper operation should also be performed after assembly.

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4.3.2 Perform the following visual inspections:

WARNING: Prior to every use the units, instruments, and accessories should be checked for damage, soiling, and completeness. All connecting cords must be checked for damage. Any inscriptions must be complete and easily legible. These inspections must be performed prior to every use.

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


4.3.3 Construction of the holding arm system

As the system can only be extended over a limited range (max. 500 mm), it is important to ensure during its installation on the operating table that all the important positions in the required working area can be reached with the holding arm. In principle, any point within the holding arm radius can be part of the working area. If possible, the arm should be positioned in advance in such a way that when released it is not necessary for the operator to catch the entire weight of the lower holding arm. The size of the operator and patient must also be taken into account here.

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4.3.4 Construction of the stand

Photo 5: First, fasten the base ① 28172 HG at the required position on a standard device rail on the operating table.

WARNING: Photo 6: All the elements of the holding system which can be tightened (e.g., levers, screws), should only be tightened manually (do not use any tools). Otherwise the elements of the system may be damaged, which could in turn endanger the patient or make it impossible to loosen the elements.

Check that the base is fixed in place correctly.

Photo 7: Insert the stand bar ② into the base ① to the required position and fix the bar in the base by tightening the fastening screw ③.

WARNING: No photo: Changing the position of the base ① and the stand bar ② also changes the position of the instrument/endoscope and/or endoscopic accessory. This could result in injury to the patient! When positioning the holding arm above the patient, the fastening screws ③ on the base (and the union nut ④ on the stand [see next page]) must not be loosened!

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
Installation and operating instructions

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4.3.5 Assembling the holding arm

Photo 8: Place the end of the stand bend (2) of the holding arm (4) over the crown teeth of the stand bar (2) and insert the stand bend into the stand as far as it will go.

Turn the stand bend (2) slightly to ensure the deepest position and correct lock-in position.

NOTE: The holding arm is delivered in one piece. It starts on one side with the stand bend, which is inserted into the stand, and ends on the other side with the handle, which is used to accept the 'KSLock' quick-release coupling, which in turn accepts the clamping jaw (11).

NOTE: The holding arm (4) and stand (2) together have the article number 28272 ECS.

Photo 9: Tighten the union nut (3) firmly. Check the correct assembly. The lock-in position is secure if it is no longer possible to turn the stand bend (2) on the retaining rod (2) when the union nut is tightened.

WARNING: First, insert the stand bend of the holding arm (2) into the stand as far as it will go. By turning the stand bend (2), the lock-in points can be clearly felt. Turn the stand bend (2) until you reach a lock-in point. Now screw the stand bend to the stand bar by screwing the union nut (3) tight manually to the right (clockwise).

While doing so, please observe the following: If the crown teeth are not correctly aligned opposite each other and the union nut (3) is tightened, this can result in the holding arm (bend) clicking in downwards when it is moved, which would loosen the union nut again and could cause an undesirable vertical position change.

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WARNING: Never perform the following procedure while the holding arm system is above the patient!

No photo: Before releasing the fastening screw ② of the base or the union nut ③ in order to reposition the holding arm system, you must support the holding arm's center of gravity so that once the fastening screws have been loosened, you can take the whole weight of the holding arm in your hand. We recommend that two people perform this procedure.

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4.3.6 Connecting the holding arm to the holding arm control unit (13)

Photo 10: Release the connector (15) from the emergency release plug (10) by retracting the cuff in the center by approx. 3 mm in the direction indicated by the arrow.

Photo 11: Then pull the connector (15) out of the emergency release plug (10).

Photo 12: Insert the connector (15) into the connection socket (14) of the holding arm control unit (11).

Photo 13: Switch the power switch (22) on the rear of the holding arm control unit to position I. The symbol (24) lights up on the front panel. The unit is now in standby mode.

Photo 14: Activate the holding arm control unit by briefly pressing the ON/OFF rocker switch (17) on the front.

The standby mode signal display (26) which was backlit until now goes out and the backlit signal display (23) indicating that the unit is ready for use now lights up, together with the holding arm display (25).

The holding arm is locked.

NOTE: If the holding arm system can be moved and/or an alarm signal is given, you should proceed according to the section "Troubleshooting and remedies" (see page 28, 29).

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WARNING: Do not bend the cable too hard. Store the cable in loops at least 15 cm in size. Cables stored on the floor (e.g., the power supply of the holding arm) must not be stepped on, squashed, or rolled over with objects. Risk of cable breaking!

WARNING: Allow the cable to trail freely. Do not wind it around the holding arm.

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4.3.7 Checking the holding force

To check the holding force, place the arm in a position used during the application.

Photo 15: Release the holding arm by pushing the two release buttons ⑥.

Then move the holding arm to the position where the spring balance is to be affixed. Then release the release buttons so that the holding arm locks.

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Photo 16: Check the holding force of the holding arm ④. To do so, affix the spring balance 28272 ECW to the handle and pull the spring balance in different directions (horizontal, vertical, diagonal). The holding force of the holding arm must be at least 20 N.

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NOTE: A force of 20 Newtons is equal to a weight of approx. 2 kg. This means that the instruments connected to the handle (e.g., endoscope, camera head) can be held securely by the holding arm system.

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Readiness for use

The EndoCrane™ holding arm system is now ready for use and can be positioned freely by simultaneously pressing the two release buttons ⑥ (described below).

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4.3.8 Functions of the handle

The release buttons ⑥ of the holding arm ④ are located on the handle ⑤. If these buttons are pressed, the holding arm is released for as long as they remain depressed and can be infinitely adjusted. Once the release buttons are released the holding arm is locked again.

The handle ⑤ has a fixation ring on its distal end (= nearest to patient), which allows free positioning of the KSLock via a ball joint in the released state (a 'Lock' imprint on the handle shows in which position the fixation ring is locked).

The 'KSLock' quick-release coupling ⑪ allows the different clamping jaws to be attached according to the application.

In our example, the clamping jaw 28272 UGN is being connected to the KSLock (see below). Alternatives are the clamping jaws 28272 UG / UK / UL / UU).

NOTE: The clamping jaws are supplied as standard. However, the connection of the clamping jaw 28272 UGN is explained for comprehension purposes.

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21

Photo 17: Press both release buttons (5) simultaneously. Then position the holding arm (4).

Then release the release buttons (5).

Photo 18: Check that the holding arm system is locked (i.e., shake it gently to see that the holding arm system cannot be moved).

Repeat the check several times.

Photo 19: The KSLock (11) is fixed firmly to the handle (5) with the ball joint.

The KSLock is a quick-release coupling system patented by KARL STORZ for the acceptance of clamping jaws.

4.3.9 Sterile covering of the holding arm

Photo 20: Preparation: Release the fixation ring and then turn the KSLock (11) until it is aligned on a straight axis with the handle (see arrow). Then turn the fixation ring in the „LOCK“ direction (the LOCK position is imprinted on the handle).

Photo 21: Pull the sterile drape over the entire holding arm system.

I.e., the drape must reach from the KSLock on one side right to the end of the stand on the other side.

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4.3.10 Fastening the clamping jaw

NOTE: The affixing of the clamping jaw is only described to aid comprehension. The clamping jaws are not supplied as standard with the holding arm system.

NOTE: In our example, we are affixing the clamping jaw 28172 UGN. Alternatives would be the clamping jaws 28272 UG / UK / UU.

Photo 23: Pierce the drape and connect the pin with the KS Lock. Simultaneous turning and pressing encourages the axis of the clamping jaw to click into place.

Photo 24: Tighten the thumb screw on the KS Lock in the 'Close' direction manually (do not use any tools). Check the correct positioning by trying to pull the clamping jaw out of the KS Lock. This should not be possible.

Photo 25: Insert the instrument into the clamping jaw (ii).

WARNING: If the clamping jaw needs to be changed during the operation, the drape must be changed before the new clamping jaw is inserted as the removal of the clamping jaw creates an unsterile opening.

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4.4 Disassembly of the EndoCrane™ holding arm system

For disassembly, proceed in the reverse order of assembly.

4.5 Important notes on avoiding damage

CAUTION: If the holding arm system is not required, it should be decommissioned. While doing so, you should aim to make the holding arm system movable in order to avoid possible damage caused by the (not permitted!) movement of lockable elements.

Always proceed as follows:

- Photo 27:** First, the holding arm should be brought unloaded into a practical, hanging position by pressing the two release buttons.
- Photo 28:** Switch off the ON/OFF rocker switch (17).
Then press and hold the ON/OFF rocker switch (17) for longer than 1 second. The operating symbol which was previously backlit and the holding arm symbol go out. The standby symbol is backlit.
- Photo 29:** Remove connector (15) from the connection socket (14).
- Photo 30:** Insert the connector (15) into the emergency release plug (16).
- Photo 31:** Press the two release buttons (6).

The holding arm system is now permanently in the released state; i.e., it can be moved and can thus be stored in the case.

Background: The connector **must** be inserted into the emergency release plug to create a movable system. To protect against mechanical damage, the emergency release plug should remain inserted during transport and storage.

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4.6 In case of fault

WARNING: In case of system failures: the emergency lock will last for at least 10 minutes! After that, the holding arm may drop! (See arrow)
Decommission the EndoCrane™ holding arm system immediately if system failures occur and remove it from the operating field.

If a fault occurs which disrupts the power supply to the holding arm while the arm is locked, the system remains locked. This emergency lock prevents the sudden dropping of the holding arm. The holding arm remains locked for at least 10 minutes with at least 15N holding force.

WARNING: If the holding arm is already in the released state (release buttons ⑥ are pressed) when the power fails or a defect affecting the power supply arises, **no** locking is triggered.
The same thing happens when the power switch is pressed and deactivated and the two release buttons ⑥ are pressed at the same time. It should be noted that the holding arm in the released state cannot be let go of and that the instrument must be removed from the site of the operation in a controlled manner.

In this case, the system must be decommissioned as described on page 27 'Procedure in case of fault'.

If the holding arm cannot be released as a result of a fault, the instruments must be disconnected from the holding arm within the time period specified above. The emergency lock provides you with a sufficient period of time to decommission the holding arm system safely in the case of a fault within 10 minutes. Once the emergency lock time expires, the holding arm may drop.

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


35

4.6.1 Consequence of power failure

The control unit emits an acoustic alarm. This is repeated at an interval of seconds.

Photo 32: All of the operating displays (23 and 24) go out. A flashing alarm display (27) is activated.

Photo 33: If the power supply is reconnected, the acoustic alarm stops and the optical alarm display (27) goes out. The operating display (26) lights up.

Photo 34: To commission the unit, the ON/OFF rocker switch (17) should be pressed briefly; the indicator for the holding arm (28) (see photo 35) lights up.

Photo 35: The alarm light (27) is reset; the operating light (26) and holding arm light (28) are on.

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
4.7 Procedure in case of fault

To decommission the holding arm, perform the following steps:

1. **Photo 36:** Remove the surgical instrument or endoscope and accessories from the clamping jaw (11).
2. **Photo 37:** Remove connector (15) from the control unit (13).
3. **Photo 38:** Connect emergency release plug (19) and connector (15).
4. **Photo 39:** Press release button (8) (holding arm (4) is released [e.g., the lock is unlocked], and the holding arm remains released).
5. **Photo 40:** Swing the holding arm (4) away from the operating table.
6. **No photo:** Disassemble the holding arm (4).

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Alarm messages, troubleshooting and remedies

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5. Alarm messages, troubleshooting and remedies

5.1 Alarm messages

The EndoCrane™ holding arm system is constantly monitored by the holding arm control unit via an integrated error detection function. A defective status is identified acoustically by an alarm signal and visually by a flashing LED in the ON/OFF rocker switch (17) field. If the alarm status is triggered, the arm must be decommissioned immediately. When troubleshooting, proceed as follows.

5.2 Troubleshooting and remedies

If any faults arise in connection with the system which cannot be resolved as described in the following table, please contact the responsible representative.

WARNING: If the faults cannot be resolved using the following table, please contact Service or send the unit in for repair. Do not attempt to repair the unit yourself!

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NOTE: Please specify the type number and serial number provided on the manufacturer's identification plate with all enquiries and in all written correspondence. On request, further documents are also available from the manufacturer.

5.3 Storage when normal operation is interrupted

If the EndoCrane™ holding arm system displays any errors which cannot be resolved using the measures specified on p. 25 – 28, the holding arm control unit must be turned off and disconnected from the power supply.

WARNING: Connecting the emergency release plug (10) to the connector (15) and then pressing the two release buttons (7) / (8) / (9) completely releases the joint lock (7) / (8) / (9)! When the holding arm is uncoupled to avoid damage and for better handling, always connect the emergency release plug to the connector!

The holding arm must be disconnected from the holding arm control unit and the holding arm must be released (i.e., the lock must be released).

To do so, the connector (15) must be pulled out of the connection socket of the control unit (photo 41) and connected to the emergency release plug (10) (photo 42), then the two release buttons (7) must be pressed (photo 43).

The holding arm (4) can now be removed from the operating table (photo 44).

The holding arm may only be checked by an authorized specialist. The system can nevertheless be stored or sent in, taking the outlined restrictions into account.

5.4 Replacing spare parts

No spare parts are available for the holding arm as repairs to the system may only be performed by KARL STORZ.

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**6. Cleaning, disinfection,
sterilization**

**6.1 Preparation of the holding
arm control unit ⑬ by the
user**

WARNING: Danger of an electric shock: avoid allowing fluids to enter the control unit's ⑬ housing. Do not store liquids on or above the device.

WARNING: Prior to all work on the control unit ⑬ (e.g., cleaning, servicing), the unit must be switched off and disconnected from the power supply.

CAUTION: No cleaning agents, abrasives, or solvents may be used in the care of the unit (control unit ⑬). Observe the manufacturer's specifications.

Disinfect the holding arm control unit ⑬ using only the 'wipe-down procedure'.

**6.2 Preparation of the holding
arm system by the user**

6.2.1 Removal of the sterile drape

WARNING: As the sterile drape for the holding arm is usually a disposable item, it must be disposed of in accordance with the applicable regulations for the disposal of medical disposables. Further information can be found in the current specialist literature. Observe the drape manufacturer's specifications.

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6.2.2 Disinfection of the holding arm and stand

Immediately after removing the holding arm from the operating table, the holding arm ④, the stand ② and the base ① should be cleaned of severe soiling with a soft, low-lint disposable cloth moistened with surface disinfectant based on quaternary ammonium compounds (surface disinfectants) or ammonium chloride (alcohol-free rapid action disinfectants) in order to rule out drying out of the contaminants.

The exterior surfaces of the unit should also be wiped clean with a low-lint disposable cloth moistened with disinfectant (same disinfectant as specified above). The disinfectant manufacturer's specifications must be observed!

WARNING: Danger of an electric shock: avoid allowing fluids to enter the control unit's housing. Do not store liquids on or above the unit.

6.2.3 Steam sterilization (not for the holding arm!)

CAUTION: Do not prepare the holding arm in steam and/or hot air sterilizers. Prior to use in the OR, cover the disinfected holding arm with a sterile drape. Only the clamping jaw and the base supplied for fastening to the operating table may be treated with the steam sterilization procedure.

The EndoCrane™ holding arm is not a sterilizable system and, as such, both sterilization using the low-temperature procedure and saturated steam must always be avoided.

Only the delivered base 28172 HG ① and the clamping jaws (28272 UL / UG / UK / UU) which are not supplied as standard should be treated using the steam sterilization procedure.

We recommend greasing the clamping jaws (with grease Art. no. 27657), in order to counteract abrasion-corrosion.

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**6.2.4 Preparation of the clamping jaws
and the base**

The clamping jaws (28272 UL / UG / UK / UU), which are not supplied as standard, are autoclavable at 134 °C +3 °C (ISO 11134), as is the base 28172 HG ①.

The following must be taken into account for the cleaning and steam sterilization.

All setting elements (knurled screws, cam levers, etc.) must be loosened in order to prevent damage via thermal expansion and to separate the surfaces for cleaning/disinfection!

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7. Service and repair

The EndoCrane™ holding arm system is maintenance-free. Prior to every use, it must be checked that the holding force is approx. 20 N in every direction. If the holding force is not achieved, the responsible Service representative should be contacted.

If any faults arise in conjunction with the system which cannot be resolved as described in the section 'Troubleshooting and remedies', please contact the responsible representative.

NOTE: Please specify the type number and serial number provided on the manufacturer's identification plate with all inquiries and in all written correspondence. On request, further documents are also available from the manufacturer.

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7.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a **year** by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

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7.2 Disposal

The control unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

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Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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In Germany you can refer repairs direct to
KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries please contact your local KARL STORZ subsidiary or authorized dealer.

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Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet the applicable laws and standards, and
- the unit has been used in accordance with its operating instructions at all times.

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**8.1 Technical data for the holding
arm control unit**

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**8.2 Interfaces to the holding arm
control unit**

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8.3 Technical data for the holding arm with stand

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8.4 Operating, storage, and transport conditions

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CAUTION: Always release the holding arm prior to shipping in order to avoid damage and for better handling (release lock). The joint locks can be completely released by connecting the emergency release plug to the connector and then pressing the release button! (See photos page 29).

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WARNING: Always lay the system's connection cord in such a way that it cannot be damaged by the personnel or instrument set.

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**8.5 Standard compliance
(for 20 7800 20)**

According to DIN 60601-1

- Type of protection against electric shock:
Protection Class I
- Degree of protection against electric shock:
Applied part type BF 

Acc. to EN 60601 / IEC 60601-1-2

Please read the Electromagnetic Compatibility Information in the Appendix (p. 41-53).

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**8.6 Directive compliance (for
20 7800 20)**

In accordance with the Medical Device Directive (MDD):

Medical device Class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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8.7 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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8.8 Software user rights

1. For KARL STORZ products with software components, the customer purchases a non-exclusive and, in the framework of these regulations, limited user right for the programmed logic, the computer programs and/or the software supplied by KARL STORZ and included in KARL STORZ products "Software".
2. The user right is only valid for the delivered form and only for use in accordance with the printed instruction manual from KARL STORZ for the products sold to the customer. Use of the software in connection with other products bought by the Customer from KARL STORZ or third parties is not permitted.
3. A transfer of the user right by the customer is only permitted if the customer informs KARL STORZ of the purchase and it is subject to the regulations included in this contract concerning the use of the software. In the case of a transfer, the customer's user right shall expire.
4. All copyrights and all rights to modify, improve, or update the software are reserved by KARL STORZ. In particular, the customer is not permitted to reproduce, copy, reverse engineer, decompile, or disassemble the software or modify it in any way. The customer is not permitted to create derivations of the software or to copy, decode, or recreate the source code of the software. § 69g Para 2 German Copyright Law (Urheberrechtsgesetz) remains unaffected.
5. Any modification or improvement of the software performed by the customer contrary to the contract shall constitute a breach of the license agreement; KARL STORZ immediately acquires all property rights and copyrights to such modifications or improvements. As a precaution, the Customer hereby grants KARL STORZ exclusive, free, and unlimited user rights to such modifications or improvements.
6. KARL STORZ recommends not using the software in connection with other software not delivered by KARL STORZ with the respective product. KARL STORZ accepts no liability for software malfunctions occurring as a result of the use of software not delivered or recommended by KARL STORZ in connection with the products delivered by KARL STORZ.
7. In the event of a failure of the Customer, or its agents, employees or representatives, to comply with the terms and conditions of the License Agreement herein granted, the user rights will terminate immediately without a specific notice of termination or any further action required by KARL STORZ or any other party.

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The holding arm system complies with EN 60601-1-2 / IEC 601-1-2 and thus fulfills the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment.

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- reorient or relocate the receiving device;
- increase the separation between the equipment;
- connect the devices to different electrical circuits.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The holding arm control unit 20 7800 20 should not be stacked directly adjacent to or on top of other devices. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The holding arm control unit 20 7800 20 is only intended for use by physicians and medical assistants. The unit complies with EN/IEC 60601 (RF emissions) and is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes. In a domestic environment the holding arm control unit 20 7800 20 may cause EMC interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the unit or filtering the connection to the public mains network.

WARNING: The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

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Table 202 Guidelines and manufacturer's declaration – electromagnetic immunity			
The holding arm control unit 20780020 and the holding arm 28272 ECS are intended for operation in an environment as described below. The user should ensure that the holding arm control unit 20780020 and the holding arm 28272 ECS are used in such an environment.			
Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	Complies ± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood or concrete or provided with ceramic tiles. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Rapid transient electric interferences/bursts acc. to IEC 61000-4-4	± 2 kV for mains lines ± 1 kV for input and output lines	Complies ± 2 kV for mains lines ± 1 kV for input and output lines	The quality of the supply voltage should be that of a typical business or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surges acc. to IEC 61000-4-5	± 1 kV push-pull voltage ± 2 kV common-mode voltage	Complies ± 1 kV push-pull voltage ± 2 kV common-mode voltage	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short-time interruptions, and fluctuations of the supply voltage acc. to IEC 61000-4-11	< 5% U_T * (> 95% dip of the U_T) for ½ cycle 40% U_T (60% dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles < 5% U_T (> 95% dip of the U_T) for 5 seconds	Conforms < 5% U_T * (> 95% dip of the U_T) for ½ cycle Conforms 40% U_T (60% dip of the U_T) for 5 cycles Conforms 70% U_T (30% dip of the U_T) for 25 cycles Conforms <5% U_T (> 95% dip of the U_T) for 5 seconds	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of the device requests continued functioning even when interruptions of the power supply occur, it is recommended to supply the device from a power supply that is free of interruptions.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	Conforms 3 A/m	The magnetic fields at the mains frequency should be at the values typical of business or hospital environments.
* Note: U_T is the mains alternating voltage before applying the test levels.			

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Table 204			
Guidelines and manufacturer's declaration – electromagnetic interference immunity for medical electrical devices that are not life-supporting			
The holding arm control unit 20780020 and the holding arm 28272 ECS are intended for operation in an environment as described below. The user should ensure that the holding arm control unit 20780020 and the holding arm 28272 ECS are used in such an environment.			
Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment - guidelines
Conducted RF disturbance acc. to IEC 61000-4-6	3 Vrms 150 kHz to < 80 MHz	3 V _{eff}	Portable and mobile RF communications equipment should be used no closer to any part of the holding arm control unit 20780020 and the holding arm 28272 ECS, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3,5/3]\sqrt{P}$ 150 kHz to < 80 MHz $d = [3,5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the nominal power of the transmitter in Watt [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. The field strength of stationary transmitters at all frequencies on site ^a should be, according to a study, less than the conformity level. ^b Interferences may occur in the vicinity of devices that carry the following symbol: 
Radiated RF disturbance acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.			
a The field strength of stationary transmitters, such as base stations of radio phones and land mobile radio services, CB radio stations, AM and FM radio, and TV stations cannot be exactly predetermined by theory. For the determination of the electromagnetic environment of the stationary transmitters, a site survey should be considered. If the measured field strength in the location at which the device is used exceeds the above conformity levels, the device should be observed to verify normal operation. If unusual performance characteristics are observed, additional measures may be necessary, such as a modified alignment or a different location for the device. b Above the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.			

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
EndoCrane™ holding arm system 28272 EC**

The holding arm control unit 20780020 and the holding arm 28272 ECS are intended for operation in an electromagnetic environment in which the interferences are controlled. The user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device, according to the output power of the communications equipment as recommended below.

Nominal power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3,5}{3} \right] \sqrt{P}$	$d = \left[\frac{3,5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters whose maximum nominal power is not listed in the above table, the recommended safety distance d in meters [m] can be determined using the equation belonging to the relevant column, where P is the maximum nominal power of the transmitter in Watt [W] as specified by the manufacturer of the transmitter.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

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WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In so doing we cannot bear the cost of transportation or the risks associated with shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:
KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

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KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Phoenicia Street, Mina Elhosh
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax: +961 1 365151

KARL STORZ Endoscopy (South Africa) (Pty) Ltd
P.O. 6061
Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Korgalynskoe shosse, 2A
Hotel „Duman“
Astana 010000, Kazakhstan
Phone/Fax: +7 7172 79 15 62, 79 15 63
E-Mail: info@karlstorz.kz

KARL STORZ Endoscopy - Gulf & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O. Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.ksegne@karlstorz-rce.com

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianle Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
06-07 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 65325548, Fax: +65 65323832
E-Mail: infoSGP@karlstorz.com.sg
E-Mail: servicesGP@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beism Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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- ① Internal part (with flexible end)
- ② Hinge
- ③ Outer sheath
- ④ Unlocking device for insert
- ⑤ Handle
- ⑥ Wheel
- ⑦ Unlocking device for internal part
- ⑧ Fixation pin
- ⑨ Fixation hole
- ⑩ Grasping bridge
- ⑪ Insert (components ①+③)

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Normal use

The CUSCHIERI Retractor is used in endoscopic surgery for the atraumatic retraction of the liver and atraumatic manipulation of other anatomical structures.

No **contraindications** relating directly to the product are currently known. Use of the instrument is contraindicated if, in the opinion of an experienced physician, the health of the patient is endangered through its use, e.g. due to the general state of health of the patient, or if the endoscopic method itself is contraindicated.

Any unauthorized change or reconstruction of the instrument is not permitted for safety reasons.

Warning: The flawless state of the instrument must be ensured before every use. If there are any signs of damage, the instrument must not be used further under any circumstances.

Warning: This instrument is **not supplied sterile** and for this reason has to be cleaned, disinfected and sterilized before its first and any other usage.

Warning: During use, the tip of the instrument must always remain within the field of view of the instrument operator.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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

Assembly

Check all the parts for damage and replace if necessary.

Preassembly of insert

Push the outer sheath (3) over the internal part (with flexible end) (1) from proximal and lock the bayonet mount with a ¼ turn.

Preassembly of handle

Fit the wheel (6) into the handle (5) by turning it until it stops, and then open by approx. 1 turn.

Final assembly

Insert the insert (11) into the handle (5) until the unlocking device for the insert engages (4).

Then, angle the hinge (2) cautiously until the unlocking device (for internal part) engages.

Turn the wheel (6) until it stops.

Check instrument for correct functioning.

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Use on the patient

Activation


Introduce the internal part (with flexible end) (1) of the assembled retractor through a suitable trocar. Using an appropriate instrument, hold the device at the grasping bridge (10) and form a 'loop' out of the flexible end. Insert the fixation pin (8) into the fixation hole (9).

The hinge (2) can also be angled if necessary using the wheel (6).

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

Deactivation

Using an appropriate instrument, hold the device at the grasping bridge (10) and pull the fixation pin (8) out of the fixation hole (9). Straighten out the flexible end of the retractor and remove it from the trocar.

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

Dismantling for cleaning

Slightly angle the hinge (2) using the wheel (6). Press the unlocking device (for internal part) (7) and stretch out the hinge (2) (do not angle any more, otherwise the lock will reengage).

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Then press the unlocking device for insert (4) and pull the insert (1) out to distal.

Turn the outer sheath (3) a ¼ turn clockwise and pull it off the internal part (with flexible end) (1).

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Cleaning, sterilization

Note: The instructions for 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail.


Caution: When preparing and using the solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration and exposure times. Prolonged exposure may result in corrosion or material damage.

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

Manual preparation

An additional manual cleaning is recommended where necessary. Clean soiling with a disinfected sponge or a disinfected soft brush. Clean the outer sheath c in an ultrasound bath and brush the inside. Following this, rinse out thoroughly with completely demineralized water and dry with sterile compressed air. Oil all movable parts at their contact points with instrument oil (Art. no. 27656) in a well-directed manner.

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Machine preparation

Ensure correct, safe positioning of the device when employing machine preparation in a washer and disinfectant (RDA). We recommend machine preparation in an RDA at up to 93 °C (+3 °C).

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Sterilization

Sterilize the instrument at 134 °C using the fractionated prevacuum method. The fully assembled instrument can be sterilized using the fractionated vacuum procedure. With other procedures, the instrument must first be dismantled. The instructions for use provided by the sterilization device manufacturer must be followed exactly.

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Operating/storage conditions

	temperature	rel. humidity
storage	-20°C ... +60°C	10% ... 90%
operation	+10°C ... +40°C	30% ... 70%

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Directive compliance

This medical device has a CE mark in accordance with Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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Normal use	5
Overview of the retractors	7
Retractors for face lifts	8
Illuminated retractors	8
50200 FR, 50200 FF	8
Manual preparation	9
Optical retractors 50200L, 50200LM, 50200LL	10
Manual preparation	10
Optical retractors with optional illumination	11
50200 R	11
Manual preparation	11
Optical retractors (DCI®) with optional illumination	12
50215 D, 50215 E	12
Manual preparation	13
Retractors for endoscopic mammaplasty and abdominoplasty	14
Standard retractors	14
50251 MR	14
Manual preparation	14
Illuminated retractors	15
496 H, 50251 R, 50251 RG	15
Manual preparation	16
Optical retractors 50200LG/LH, Man. preparation	17
Optical retractors with optional illumination	18
50251 LD, 50251 LC	18
Manual preparation	18
Optical retractors (DCI®) with optional illumination	19
50216 D	19
Manual preparation	20
Retractors for reconstructive surgery	21
Standard retractors	21
50205 EK	21
Manual preparation	21
Illuminated retractors	22
496 H, 50251 R, 50251 RG	22
Manual preparation	22
Optical retractors	23
50251 LG, 50251 LH	23
Manual preparation	23
Optical retractors with optional illumination	24
50200 L, 50200 LM, 50200 LC, 50200 LD	24
Manual preparation	24
Cleaning, disinfection, sterilization	25
Subsidiaries	27

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Intended use

The retractors are used in plastic surgery for creating and maintaining working space required for the operation in question.

Depending on the model used, endoscopes can also be used for optical visualization of the operating field and fiberoptic light carriers can be employed for its illumination.

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No **contraindications** relating directly to the product are currently known. The use of the product is contraindicated if, in the opinion of the responsible physician, such an application could cause danger to the patient, e.g., due to the general condition of the patient, or if the endoscopic method as such is contraindicated.

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Application

WARNING: A risk of injury to the patient/user and/or damage to the products exists if the instruction manuals and warning symbols on the products are not observed. Observe the instruction manual and warnings.

WARNING: Risk of injury to the patient/user and/or damage to the products: The instruction manuals and interface specifications for the medical devices used in combination must be observed precisely.

WARNING: Risk of injury to the patient/user: Overloading the instrument (as a result of applying or exerting too much force on it), may cause the instrument to break and/or bend and cause functional/operational malfunctions, which can result in the injury of the patient or user. Never overload the instrument. Do not bend bent instruments back to their original state (tensile fractures cause corrosion).

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
WARNING: Risk of infection: These instruments are **not** sterile when delivered and must therefore be cleaned, disinfected and/or sterilized prior to initial use and each subsequent reuse.

WARNING: Risk of injury to the patient/user and/or damage to the products: Combinations of medical devices are only assured to be safe if they are identified as such in the respective instruction manuals or the intended purpose and interface specifications of the devices used in combination permit this.

WARNING: Risk of injury to the patient: Never rest the end of a light cable on or under a sterile drape. The intensity of the light can cause burns to the patient and/or the sterile drape.

WARNING: Risk of injury to the patient: Lighting power should be kept as low as possible and used for as short periods of time as possible to avoid potential tissue damage due to heating.

WARNING: Risk of injury to the patient/user: Ensure that the instruments are in perfect working order before each use. If there are any signs of damage, an instrument must not be used under any circumstances.

WARNING: Follow the application and storage conditions in this manual. Only use channels for the designed purpose.

WARNING: Risk of injury to the user: The retractor should be handled carefully during its application and preparation to avoid injuries.

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Manual preparation 'illuminated retractors'

Remove light carriers from the retractor by pressing the corresponding slide.

50200 FR (1st photo): Push the slide gently in the direction indicated by the arrow using your thumb.

50200 FF (2nd photo): Place the thumb and index finger as indicated in the image on the left. Then, using your thumb, press gently from below in the direction indicated by the arrow so that the light carrier is freed.

CAUTION: The index finger stops the light carrier from sliding out of the retractor too quickly and thus avoids the possibility of it falling on the floor.

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Parts without an integrated light cable may be cleaned in an ultrasound bath (maximum five minutes). Finally, rinse with water so that all cleaning residue is removed.

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Immerse all parts completely in a disinfectant solution making sure that all air bubbles are removed and taking into account the instructions of the chemical manufacturer concerning concentration and exposure times. Rinse thoroughly with microbiologically pure/sterile water and dry with a lint-free disposable cloth or sterile compressed air.

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
**Optical retractors with
optional illumination**

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50200 R: Optical retractor, with handpiece, width of distal end 30 mm, working length 10 cm, with fixing screw for regulation of the insertion depth of the endoscope, for use with **HOPKINS®** telescopes of 4 mm diameter or with the fiberoptic light carrier 50200 RL for non-endoscopic applications.

- ⑤ Handpiece
- ⑥ 50230 BA **HOPKINS®** telescope 30°, diameter 4 mm
- ⑦ Fiberoptic light carrier

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Manual preparation 'Standards retractors'

Clean the surfaces carefully with a disinfected sponge or a disinfected soft brush. Brush any lumina well and/or rinse through.

Parts without an integrated light cable may be cleaned in an ultrasound bath (maximum five minutes). Finally, rinse with water, so that all cleaning residue is removed.

Immerse all the parts in a disinfectant solution making sure that all air bubbles have been removed and that the parts are completely covered, and taking into account the instructions provided by the chemical manufacturer on concentration and exposure times. Rinse thoroughly with microbiologically pure/sterile water and dry with a lint-free disposable cloth or sterile compressed air.

NOTE: The procedures for machine cleaning and sterilization are described on page 25.

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
**Retractors for
plastic surgery**

 Retractors for endoscopic mammaplasty
and abdominoplasty

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
496 H: Illuminated retractor with fiberoptic light carrier, blade width 25 mm, length 14 cm

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- ⑬ Atraumatic rounding
- ⑭ Light carrier

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
50251 R: Illuminated retractor with fiberoptic light carrier, blade width 30 mm, length 9 cm, with prongs, with smoke evacuation channel

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- ⑮ Prongs for a secure hold
- ⑯ Light carrier
- ⑰ Smoke evacuation channel

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50251 RG: Illuminated retractor with fiberoptic light carrier, blade width 35 mm, length 12 cm, with atraumatic prongs, with smoke evacuation channel

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- ⑱ Atraumatic prongs
- ⑲ Light carrier
- ⑳ Smoke evacuation channel

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**Manual preparation 'Illuminated
retractors'**

Remove light carrier from the retractor. This is done by sliding the light carrier out of the guide laterally.

Clean the surfaces carefully with a disinfected sponge or a disinfected soft brush. Brush any lumina well and/or rinse through.

CAUTION: Parts with an integrated light cable must not be cleaned using ultrasound, as this may damage them.

Parts without an integrated light cable may be cleaned in an ultrasound bath (maximum five minutes). Finally, rinse with water, so that all cleaning residue is removed.

Immerse all the parts in a disinfectant solution making sure that all air bubbles have been removed and that the parts are completely covered, and taking into account the instructions provided by the chemical manufacturer on concentration and exposure times. Rinse thoroughly with microbiologically pure/sterile water and dry with a lint-free disposable cloth or sterile compressed air.

NOTE: The procedures for machine cleaning and sterilization are described on page 25.

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Retractors for plastic surgery

Retractors for endoscopic mammoplasty
and abdominoplasty

Optical retractors

50251 LG: Optical retractor, for mammoplasty and abdominoplasty, blade width 25 mm, with two individually settable lateral smoke evacuation channels, for use with **HOPKINS®** telescope 50253 BA

① Smoke evacuation channels
② 50253 BA **HOPKINS®** telescope 30°, diameter 10 mm

50251 LH: Optical retractor, for breast and reconstructive surgery, for creating an operative pocket, with handpiece, laterally oblique, with lateral telescope guide, for use with **HOPKINS®** telescope 50230 BA

③ Ergonomic design
④ 50230 BA **HOPKINS®** telescope 30°, diameter 4 mm

Manual preparation 'Optical retractors'

Remove the telescope from the retractor.

Clean the surfaces carefully with a disinfected sponge or a disinfected soft brush. Brush any lumina well and/or rinse through. Parts without an integrated light cable may be cleaned in an ultrasound bath (maximum five minutes). Finally, rinse with water, so that all cleaning residue is removed.

Immerse all the parts in a disinfectant solution making sure that all air bubbles have been removed and that the parts are completely covered, and taking into account the instructions provided by the chemical manufacturer on concentration and exposure times.

Rinse thoroughly with microbiologically pure/sterile water and dry with a lint-free disposable cloth or sterile compressed air.

① **NOTE:** The procedures for machine cleaning and sterilization are described on page 25.

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Optical retractors with optional illumination

50251 LD: Optical retractor, for creating an operative pocket, with two separate smoke evacuation channels, with handpiece for single-handed operation, for use with **HOPKINS®** telescope 50251 BA or with the fiberoptic light carrier 50251 LT for non-endoscopic applications.

50251 LC: Same, without smoke evacuation channels

- ②③ Retractor with ergonomic design
- ②④ 50251 BA **HOPKINS®** telescope
- ②⑤ Fiberoptic light carrier
- ②⑥ Smoke evacuation channels
- ②⑦ Guide sheath

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Retractors for reconstructive surgery

Standard retractors

50205 EK: Tissue retractor, with catch, for the
retraction of tissue during surgical interventions on
the face.

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Manual preparation 'Standard retractors' for reconstructive surgery

Clean the surfaces carefully with a disinfected
sponge or a disinfected soft brush. Brush any
lumina well and/or rinse through.

Parts without an integrated light cable may be
cleaned in an ultrasound bath (maximum five min-
utes). Finally, rinse with water, so that all cleaning
residue is removed.

Immerse all the parts in a disinfectant solution
making sure that all air bubbles have been re-
moved and that the parts are completely covered,
and taking into account the instructions provided
by the chemical manufacturer on concentration
and exposure times.
Rinse thoroughly with microbiologically pure/sterile
water and dry with a lint-free disposable cloth or
sterile compressed air.

NOTE: The procedures for machine
cleaning and sterilization are described on
page 25.

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
Illuminated retractors

496 H: See Section 'Retractors for endoscopic mammaplasty and abdominoplasty' page 15

50251 R: See Section 'Retractors for endoscopic mammaplasty and abdominoplasty' page 15

50251 RG: See Section 'Retractors for endoscopic mammaplasty and abdominoplasty' page 15

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Retractors with integrated illumination

50251 RBB/RSB: Retractor for the endoscopic elevation of the latissimus dorsi muscle with integrated illumination. Can be used with the HOPKINS® II enlarged view forward oblique telescope 30° 26105 BA.

50251 RHB: Handle to be used with 50241 RBB/RSB

Assembly

Connect the light cable to light port ①.
Guide handle ② over the light cable and screw it onto thread ③ on the retractor by twisting.
Telescope sheath 50200 LS 4 can be optionally pushed over the attachment into the retractor and secured by rotating one quarter.

Disassembly

To disassemble, remove the light cable from light port ①. Remove the handle by unscrewing it from retractor thread ③.
Unfasten the telescope sheath by turning one quarter and remove telescope sheath ④ from the retractor. The sheath should also be detached for cleaning purposes.

Manual (pre-)cleaning

Carefully clean the surfaces with a sponge or a soft brush. Brush any lumina through. Then rinse all the parts thoroughly with completely demineralized water and dry with sterile compressed air.

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Cleaning, disinfection and sterilization

CAUTION: The instructions for 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail. The current list of KARL STORZ approved materials can be found on the Internet at 'www.karlstorz.com'.

CAUTION: When producing and using cleaning and disinfecting solutions, the manufacturers' instructions for concentration and exposure times must be followed precisely.

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Manual preparation

The type of manual preparation employed depends on the respective instrument group. For this reason, the manual preparation stages are described in the sections above for each of the relevant instrument groups.

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Machine preparation

Before machine preparation, manually perform a preliminary cleaning (possibly with ultrasound bath).
Ensure that all parts are stored correctly to avoid any damage occurring when machine cleaning is in progress. To this end, the parts should ideally be placed in a tray holder.
The following procedures may be used:
• chemical disinfection
• thermal cleaning at 93°C

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Sterilization

These instruments can be gas sterilized or steam sterilized at 134°C + 3°C. The instrument can be sterilized via the fractionated vacuum method in its assembled state. For other sterilization methods, it is necessary to dismantle the instrument. Follow the sterilization unit and disinfectant manufacturers' instructions exactly.

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Disposal

To dispose of the retractors, no special measures are necessary. National laws/regulations must be observed.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kse.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kse.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520 México México
Phone: +52 55 5250-5607, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufralegui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtravägen 14,
12724 Skärholmens, Sweden
Postal address: P.O. Box 8013,
14108 Kungälv, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4, kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlystvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopy Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stoplep.be

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta 1
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Coloman, nr. 74, Sector 4
041393 Bucuresti, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Gouran 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

OOO KARL STORZ Endoscopy — VOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-vostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391
E-Mail: info@karlstorz-erng.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kazakhstan
Sary-Arka str., 6, BC "Arman", off. 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) — Units 205 — 206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-erng.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanrong Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg


KARL STORZ Endoscopy Korea Co. Ltd.
Jeilong Tower 1, Suite 811
17-1 Jeongja-dong
Seongnam, 463-847, Korea
Phone: +82-70-3450-7474, Fax: +82-70-8277-3299

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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
- ① B.I.O.H. outer sheath 4 mm (26252BO)
- ② Compact hysteroscope 26252BB/BH (single, short/long handle without outer sheath and valves)
- ③ Valve connector for connecting to tubing set 031317-10 (mtp) or Monobloc (26252 BM)
- ④a Monobloc adaptor (26252 BM) to connect the B.I.O.H. compact hysteroscope to a tubing set
- ④b Tubing set 031317-10 (mtp)
- ⑤ Light port
- ⑥ Working channel
- ⑦ Sealing cap (7919291)
- ⑧ Suction and irrigation valves (26252BV)

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1. Intended use

The BETTOCCHI B.I.O.H. compact hysteroscope is used in combination with endoscopic accessories for diagnosis and/or treatment in hysteroscopy.

No **contraindications** relating directly to the product are currently known. The BETTOCCHI B.I.O.H. compact hysteroscope must not be used if, in the opinion of the attending physician, such an application would endanger the patient, e.g., due to the general condition of the patient.

Unauthorized conversions or modifications to the instrument are prohibited for safety reasons.

2. Safety instructions

! WARNING: Risk of injury/material damage if the instruction manual/warnings are not observed: **Read and observe the instruction manual and warnings for the devices and the devices used in combination.**

! WARNING: Risk of infection for patient/staff: These instruments are not sterile when delivered. Clean, disinfect and sterilize the instruments before initial use and before every subsequent use using validated procedures.

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**BETTOCCHI B.I.O.H.
compact hysteroscope**

! WARNING: Risk of injury to the patient/user and/or damage to the products: Combinations of medical devices are only assured to be safe if they are identified as such in the respective instruction manuals or if the intended purpose and interface specifications of the devices used in combination permit this.

! WARNING: Risk of injury for patient: The incorrect fitting of an instrument's assembled components can result in injury to the patient. Check the secure fitting of the components before each application by pulling them gently in both directions and/or turning them without exerting excessive force.

! WARNING: Risk of injury for patient/user: Do not use damaged instruments/accessories (e.g., if they have rough surfaces, sharp corners, burred edges, or projecting parts). Always check instruments/accessories before and after use to ensure that nothing is missing. Test for proper operation. Replace incomplete instruments/accessories.

! WARNING: Risk of infection for patient/staff: Check instruments for visible contamination (blood, tissue remnants). This indicates that preparation has not been performed, or has been performed incorrectly. In such cases, clean, disinfect and sterilize the instrument again using validated procedures.

! WARNING: Risk of infection for patient/staff as a result of incorrect preparation: Follow the "Cleaning, sterilization and care of KARL STORZ instruments" instructions.

! WARNING: Risk of injury to the patient/user: The instruments may only be used in accordance with the medical rules and procedures for endoscopy.

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! WARNING: Risk of injury to the patient/user: Overloading the instrument (as a result of applying or exerting too much force on it) may cause it to break, bend and malfunction and consequently injure patients or users. Do not place too much pressure on the instrument. Do not attempt to bend bent instruments back to their original position (tensile cracks lead to corrosion).

! WARNING: Risk of injury: Only employ the instrument if you possess an appropriate medical qualification and are acquainted with the application.

! WARNING: Risk of burns for patient/user! The high level of light energy means that the distal end, the light port and adjacent components may become so hot that they can cause burns! Failing to maintain sufficient distance between the light outlet and tissue can result in tissue damage. Avoid direct tissue contact with the distal end and with the light connection of the telescope and the light cable. Never rest the endoscope or the light cable on the patient or where it is in direct contact with the operating accessories. Always select the lowest light output of the cold light source which still allows optimal illumination of the operating field. Please avoid using the light for long periods at a high intensity.

! WARNING: Risk of infection or damage to products: Effective preparation cannot be guaranteed if non-validated procedures are employed. Also refer to the list of approved agents for preparation. This can be found on the Internet (www.karlstorz.com) and in the 'Cleaning, sterilization, disinfection and care of KARL STORZ instruments' instructions.

! WARNING: Risk of infection for patient/user: The successful sterilization of a medical device can only be ensured when preparation is performed under aseptic conditions.

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 **WARNING:** Danger of glare! Never look into the free end of either a connected fiberoptic light cable or endoscope.

 **WARNING:** If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image becomes cloudy or if no image or only parts thereof are visible.

 **WARNING:** Avoid high irrigation pressures during application: Risk of pyrogen effect!

 **CAUTION:** Observe the application and storage conditions. Channels must only be used for their intended purpose.

 **NOTE:** Sealing systems must only be used once, as otherwise the impermeability can no longer be guaranteed.

 **NOTE:** If necessary, clean the lens and eye lens with a cloth soaked in alcohol.

 **NOTE:** The loss of certain functions (e.g., illumination failure during an intervention) can result indirectly in risks for the patient. For this reason, an additional, ready-to-use lamp and additional spare device should always be kept ready.

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3. Cleaning, disinfection, sterilization

 **CAUTION:** The manual 'Cleaning, sterilization and care of KARL STORZ instruments' must be followed. The cleaning, disinfection and sterilization procedures are explained therein in detail. The list of KARL STORZ approved materials can be found on the Internet at 'www.karlstorz.com' or in the manual 'Cleaning, sterilization and care of KARL STORZ instruments'.

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3.1 Preparation for cleaning

1. If necessary, remove instruments from the working channel.
2. Unscrew the light cable and light cable adaptor.
3. Remove the valves for cleaning (when not in the locked position) from the handle by pulling sharply.
4. Inspect all O-Rings (total: 6) for damage and replace if necessary.

NOTE: Make note of the three different sizes of the O-Rings for the valve 26252BV.

5. Replace damaged sealing rings.
6. Remove the outer sheath ① by holding the working insert ① on the union nut and pulling it forward.
The outer sheath ① engages with a clicking sound and can be pulled off toward the front.
7. Inspect the O-Ring from the working channel and replace if necessary.
8. Remove the sealing cap ④ from the compact hysteroscope ②, check again and replace if necessary.

The following **replacement seals** can be found in the spare parts box, art no. 26252SP:

10 x 5921100 O-Ring (for valve 26252BV)
10 x 5922700 O-Ring (for valve 26252BV)
10 x 5936400 O-Ring (for valve 26252BV)
5 x 7919291 Sealing cap (for working channel)
5 x 5922000 O-Ring (for outer sheath 26252BO)

9. Press the locking mechanism on the handle and remove the Monobloc adaptor ⑦ from the valve connector ③.
10. Check the O-Ring of the Monobloc adaptor ⑦, and replace if necessary.
11. Place the light cable adaptor, the sealing cap (7919291) and Monobloc adaptor ⑦ in the small parts tray of the tray holder for the final machine cleaning and disinfection.

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
3.2.2 Compact hysteroscope (26252BB/BH)

1. Thoroughly clean (with a disinfected brush) on the handle of the B.I.O.H. shaft all (lumen)
 - both lumina of the Monobloc connector ③ (brush art. no.: 27650 B),
 - the lumen of the working channel ④ (brush art. no.: 27651 UA),
 - both lumina of the B.I.O.H. handle which are visible when the suction and irrigation valves ⑧ are removed.

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2. Place the cleaning adaptor 39501 XCA in the suction and irrigation ports (with suction and irrigation valves removed) and thoroughly rinse the lumen of the suction and irrigation connectors using a rinsing aid (e.g., syringe).
3. Rinse through (with microbiologically pure/sterile water) on the handle of the B.I.O.H. shaft all (lumen)
 - both lumina of the Monobloc connector ③
 - the lumen of the working channel ④
 - both lumina of the B.I.O.H. handle which are visible when the suction and irrigation valves ⑧ are removed.
4. Then dry all parts (with sterile compressed air) and on the handle of the B.I.O.H. sheath all (lumen)
 - both lumina of the Monobloc connector ③
 - the lumen of the working channel ④

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
• both lumina of the B.I.O.H. handle which are visible when the suction and irrigation valves ⑧ are removed.

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
3.5 Sterilization

Steam sterilize all parts with moist heat at 134°C, hold time 5 minutes, using the fractionated vacuum method. With the fractionated vacuum method, this instrument can be sterilized in the assembled state, apart from the sealing systems (these must be dismantled before being sterilized).

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NOTE: Follow the instructions for use provided by the manufacturer of the sterilization device.

4. Assembly

1. Screw the light cable adaptor on the light inlet piece ①.
2. Click on the outer sheath.
3. Insert the valves in the handle.
4. Place the sealing cap on the working channel.
5. Check that the instrument is functioning correctly.

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	Temperature	Rel. humidity
Storage	-20°C ... 60°C	10% ... 90%
Application	10°C ... 40°C	30% ... 70%

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
STORZ
KARL STORZ—ENDOSKOPE

Take-apart instruments with sliding shaft (CLEANFIT® lug-and-recess locking mechanism)

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97000003

- ① Front of a disassembled instrument with lug-and-recess locking mechanism
- ② Back of an assembled instrument with lug-and-recess locking mechanism and locking bar
- ③ Lug
- ④ Slide
- ⑤ Body of instrument (distal with working tip)
- ⑥ Guide
- ⑦ Recess for inserting the guide (not visible in the photograph because it is located in the body of the instrument)
- ⑧ Small ring component
- ⑨ Body of instrument (proximal with ring)
- ⑩ Safety bar
- ⑪ Opening aid 459999

1 Intended use

The instruments are used for gripping, removing, cutting and dissecting skin, tissue, cartilage, bone and tamponades in ENT and in skull base and endoscopic pituitary gland surgery.

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2 Contraindications


No contraindications relating to the product are currently known. Use of the instruments is contraindicated if, in the opinion of the responsible physician, the patient's health is endangered through their use, e.g., due to the patient's general condition or if the endoscopic method as such is contraindicated.

Unauthorized conversions or modifications to the instruments are not permitted for safety reasons

3 Safety instructions

WARNING: Risk of injury: Only employ the instruments if you possess an appropriate medical qualification and are acquainted with the application.

WARNING: Risk of injury/material damage if the instruction manual/warnings are not observed: Read and observe the instruction manual and warnings for the products and the products used in combination.

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Symbol description

-  **WARNING:** Failure to observe may result in injury or even death.
-  **CAUTION:** Failure to observe may result in damage to or even destruction of the product.
-  **NOTE:** Special information on the operation of the instrument.
-  Read accompanying documents before use
-  Manufacturer

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-  **WARNING:** Risk of infection: KARL STORZ instruments are provided non-sterile and must be cleaned and sterilized/high-level disinfected before initial use and before every subsequent use.
-  **WARNING:** Risk of infection for patient: Prior to sterilization inspect the instrument for visible soiling (blood, tissue remnants), which indicates that cleaning and disinfection have not been performed or have been performed incorrectly. Repeat cleaning as necessary.
-  **WARNING:** Risk of infection for patient/staff: Ensure the sterility of the instruments and surroundings.
-  **WARNING:** Risk of infection or damage to products: Effective reprocessing cannot be guaranteed if non-validated procedures are employed.
-  **WARNING:** Risk of infection from unsterilized instruments. Only use the sterilization methods recommended in the instruction manual.
-  **WARNING:** Risk of injury to the patient/user and/or damage to the products: Combinations of medical products are only assured to be safe if they are identified as such in the respective instruction manuals or the intended purpose and interface specifications of the products used in combination permit this.
-  **WARNING:** Risk of injury to the patient/user: Damaged or incomplete instruments/accessories may cause malfunctions. Always check instruments/accessories before and after use to ensure that nothing is missing/damaged. Replace incomplete/damaged instruments/accessories and do not continue to use them. Before every use check that all of the instruments required for the operation are functioning correctly.

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


- WARNING:** The product described here must be used in accordance with the medical rules and procedures of endoscopy recognized for this method.
- WARNING:** Check the instrument and accessories used in combination before every application for unintentional rough surfaces, sharp corners, burred edges, and projecting parts, which can endanger personal safety. Damaged, bent, or incomplete instruments or accessories must not be used further under any circumstances.
- WARNING:** Risk of injury to the patient/user: Overloading the instrument (as a result of applying or exerting too much force on it) may cause it to break, bend, and malfunction, and consequently injure the patient or user. Do not exert too much pressure on the instrument. Do not attempt to bend bent instruments back to their original position (tensile cracks lead to corrosion).
- CAUTION:** Observe the application and storage specifications stated in the instructions.

4 Disassembly

- The instruments must be disassembled for reprocessing.
- If the instrument has a safety bar, open this bar on the back of the instrument in the direction indicated by the arrow. It must be possible to do this without using force (see Fig. A).
- Position the tip of your thumb on the instrument exactly as shown in Fig. B. It is very important that you do not place the tip of your thumb over the lug, as this would prevent the instrument from opening in the next step of disassembly. The tip of your thumb must be positioned on the edge next to the lug (see arrow in Fig. B).
- Press gently against the edge (see Fig. C) with the tip of your thumb and then press the slide downward.
 - The locking mechanism (the lug) is now open (see Fig. D).
- The small ring component can now also be removed (see Fig. E).

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



4.1 Disassembly using the opening aid 459999

CAUTION: Carefully apply the opening aid. To avoid bending the take-apart instrument, do not apply the opening aid 459999 using excessive force.

1. Hold the opening aid ⑪ parallel to the take-apart instrument as shown in Fig. F.
2. Push the opening aid over the instrument until the notch of the opening aid lies directly on the lug ③. To facilitate handling, open the grip of the instrument and press the opening aid gently against the edge of the upper sheath half (see Fig. G).
3. The locking mechanism (the lug) is now open (see Fig. H).
4. The small ring component ④ now also be removed (see Fig. I).

5 Assembly

1. Re-attach the disassembled ring component to the body of the instrument. The grip clicks audibly into place in the body of the instrument (see Fig. J).
2. There is a guide in the middle of the slide of the instrument and a corresponding recess in the middle of the body of the instrument. Press the guide into the recess (see arrow in the photograph, see Fig. J).
3. Press the slide of the instrument so that the circular notch is located below the lug (see arrow in the photograph, see Fig. K).
4. Press the circular notch of the slide into the lug. This may require moving the grip of the instrument to and fro a little so that the circular notch of the slide can engage in the lug. Ensure that the notch is fully engaged, as otherwise the safety bar cannot be closed in the next step of assembly.
5. Close the safety bar (see arrow direction in Fig. L).

6 Reprocessing

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

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CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing.

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6.1 Accessories

Accessories required for carrying out reprocessing

- Brushes: 27652
- Opening aid: 459999

6.2 Disassembly

The medical device must be opened as far as possible prior to cleaning and sterilization/high-level disinfection (see Chapter 4).

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6.3 Pre-cleaning at the point of use

1. Immediately after a procedure, at the point of use (e.g. operating room), wipe the instruments with a soft, lint-free, disposable cloth, moistened with distilled/demineralized water or cleaning detergent to remove all heavy soiling.
2. Contain and transport the contaminated instruments to the decontamination area for manual cleaning.
3. If immediate cleaning is not feasible, place the instruments in a container and soak in neutral/ mild pH enzymatic cleaning solution (e.g. Enzofol), diluted with distilled/demineralized water to proper concentration per manufacturer's instructions.

6.4 Manual cleaning

1. Disassemble the instruments.
2. Thoroughly rinse instruments with distilled/ demineralized water to remove gross debris.
3. Completely immerse instruments in a mild/neutral pH enzymatic solution. Keep immersed for 10 minutes.
4. While immersed, brush the instrument with the cleaning brush.
5. After immersion period, rinse the instruments with distilled/demineralized water.
6. Dry the instruments with a lint-free soft cloth or compressed air.

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6.5 Inspection

The cleaned instrument must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Joints and glide surfaces must be lubricated locally with instrument oil prior to steam sterilization.
- Afterwards, a functional check must be carried out.

NOTE: The oil used for this purpose must be suitable for the subsequent sterilization process (silicone-free and paraffin- or white oil-based).

6.6 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 -10, ISO 11607 Parts 1 + 2, DIN 58953).

6.7 Sterilization

The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following steam sterilization parameters have been validated by KARL STORZ for use in the U.S.

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7 Application and storage conditions

	Temperature	Relative humidity
Storage	-20 °C ... +60 °C	10 % ... 90 %
Application	+10 °C ... +40 °C	30 % ... 70 %

8 Directive compliance

This medical device bears a CE mark in accordance with MDD 93/42/EEC. The four-digit code number after the CE mark indicates the notified body involved in certification in Europe.

9 Disposal

For disposal, no special measures are necessary. Disposal must be carried out in accordance with the respective local and national laws and regulations.

10 Packaging symbols

For the meanings of the symbols printed on the packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. It can be downloaded from www.karlstorz.com.

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Hersteller

Symbol description

WARNING: Failure to observe may result in injury or even death.

CAUTION: Failure to observe may result in damage to or even destruction of the product.

NOTE: Special information on the operation of the instrument.

Read accompanying documents before use

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
- ▷ Loading cone
- ◁ Suction head
- △ Plunger
- ▽ Fixation screw
- ▷ Trigger
- ◁ Handle
- △ Suction line connection
- ▽ Suction inlet

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1. Intended use

Hemorrhoid ligature instruments are used to place ligature rings over a hemorrhoid in order to cut off its blood supply. The hemorrhoid ligature instruments are loaded with a ligature ring with the help of a cone.

No **contraindications** relating directly to the product are currently known. The suction ligature instrument must not be used if, in the opinion of the responsible physician, the health of the patient is endangered through such use, e.g., due to the general condition of the patient.

Unauthorized conversions or modifications to the instrument are prohibited for safety-related reasons.

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2. Safety instructions

! WARNING: Risk of injury to patient: Warnings and instruction manuals for the products being used in combination must be observed. Incorrect product combinations can result in injury to the patient and/or user, or damage to the products. Combined use only in compliance with the intended purpose.

! WARNING: Risk of injury to patient: Damaged or incomplete instruments/accessories may not function normally. Check instruments/accessories are complete and are functioning correctly before and after each use. Replace damaged/incomplete instruments/accessories and do not reuse them.

! WARNING: Risk of injury to patient: Too sharp or burred edges and corners can injure the patient. Check instruments before and after each use for sharp or burred edges. Replace damaged instruments/accessories and do not reuse them.

! WARNING: Risk of injury to patient: The incorrect fitting of an instrument's assembled components can result in injury to the patient. Before each use, the secure position of the components must be checked by gently pulling and/or turning them.

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WARNING: Risk of injury: The instruments described here must be used in accordance with the medical rules and procedures recognized for endoscopy.

WARNING: Risk of injury to patient/surgeon: Overloading the instrument (as a result of applying or exerting too much force on it), may cause the instrument to break, bend, and malfunction, and consequently injure patients or users. For this reason, never overload the instrument. Bent instruments must not be bent back into shape. Instead, these instruments must be disposed of.

WARNING: Risk of injury to patient/surgeon: The instruments may only be used by persons with an appropriate medical qualification and who are acquainted with the technique.

WARNING: Risk of infection to staff/patients: These instruments are not sterile when delivered. They must be cleaned, disinfected and sterilized before the first use and before every subsequent use using validated procedures. Ensure that the instruments and surroundings are sterile.

WARNING: Risk of infection to staff/patients: Check the instrument for any visible soiling (blood, tissue remnants), which indicates preparation was not carried out or was carried out incorrectly. If necessary, prepare again using validated procedures.

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CAUTION: The 'Cleaning, sterilization, and care of KARL STORZ instruments' instructions must be followed. A list of KARL STORZ approved materials can be found on the Internet (www.karlstorz.com) or in the manual 'Cleaning, sterilization, and care of KARL STORZ instruments'.

NOTE: Follow the application and storage conditions in this manual!

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3. Assembly

1. Connect the suction ligature instrument to your vacuum pump with a suitable tube.

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4. Performing the ligation

2. Turn the vacuum pump on by pressing the power switch.
3. Use the proctoscope to advance the suction ligature instrument towards the hemorrhoid.
4. Block off the suction inlet above the handle using your thumb.
- This creates an underpressure, which sucks the hemorrhoid into the suction head.
5. Press the trigger.
- The plunger pushes the ligature ring off the suction head and positions it over the suctioned hemorrhoid.
6. To remove the underpressure, remove your thumb from the suction inlet.

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5. Disassembly

1. Unscrew the suction head.
2. Unscrew the fixation screw and remove the release mechanism.

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6. Cleaning, disinfection, and sterilization (preparation)

 **WARNING:** Risk of infection: The methods specified in the instruction manual have been validated by KARL STORZ. Only use these methods, as otherwise sterility cannot be guaranteed.

 **WARNING:** Product damage: The preparation method must be validated by the user on site.

 **WARNING:** Product damage: The methods specified in the instruction manual have been validated by KARL STORZ. Only use these procedures.

6.1 Manual (pre)cleaning

 **NOTE:** We recommend washing and rinsing the instruments thoroughly immediately after clinical applications to avoid contaminants adhering.

 **NOTE:** They must be transported to the preparation location immediately.

1. Clean the surfaces carefully immediately after use with a moistened, low lint disposable cloth (preferably with distilled water).

6.2 Manual cleaning

 **NOTE:** We recommend performing an ultrasound treatment at 35 kHz (5 minutes).

 **NOTE:** Refer to the chemical manufacturers' instructions on concentration and exposure times.

1. Immerse all parts in a cleaning solution making sure that all air bubbles have been removed and that the parts are completely covered.

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Suction ligature instrument

2. Brush all parts with a soft, disinfected sponge/brush, taking care to avoid spraying.
3. Then flush the cleaning solution through the suction line connection using an injection aid/cleaning gun.
4. Rinse all parts thoroughly with microbiologically pure/sterile water.

6.3 Manual disinfection

NOTE: Refer to the chemical manufacturers' instructions on concentration and exposure times.

NOTE: Observe the health and safety regulations.

1. Immerse all parts in a disinfectant solution (preferably VAH/DGHM-approved) making sure that all air bubbles have been removed and that the parts are completely covered.
2. Then flush the disinfectant solution through the suction line connection using an injection aid/cleaning gun.
3. Rinse all parts thoroughly with microbiologically pure/sterile water.
4. Finally, dry all parts with a low lint disposable cloth / sterile compressed air.

6.4 Machine cleaning and disinfection

WARNING: Risk of infection: In cases of severe, optically contaminated medical devices, manual cleaning must be performed in advance.

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CAUTION: Risk of material damage: If the individual parts are not placed safely in the washer and disinfecter, this can lead to damage during machine preparation or small parts can be lost. Ensure the safe positioning of all parts and store the small parts in the small parts tray.

1. Place all parts correctly in a suitable tray holder.
2. Place the tray holder with its contents in the washer and disinfecter.
3. Connect the suction line connection of the handle up to the washer and disinfecter.
4. Prepare the instrument set in the washer and disinfecter with thermal disinfection at up to 93°C (alkaline cleaning agent ≤ 11.5 pH, please refer to the list of approved chemicals at www.karlstorz.com).

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7. Care

NOTE: Lubricate all the parts which move against each other with a drop of white paraffin oil-based care oil.

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8. Sterilization
» Performing sterilization

NOTE: Use the fractionated vacuum procedure (acc. to DIN EN 550/ISO 11134) observing all nationally applicable regulations.

WARNING: Risk of infection: Only sterilize cleaned, disinfected, dried, and functioning medical devices in suitable packaging.

NOTE: Follow the instructions for use given by the manufacturer of the sterilization device.

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Suction ligature instrument

NOTE: The instrument can be sterilized without being disassembled in the fractionated prevacuum procedure.

1. Steam sterilize the instrument at 134°C (device-specific) using the fractionated prevacuum procedure (5 minutes hold time).

9. Operating/storage conditions

CAUTION: Risk of material damage: Store the sterilized instruments in sterile goods packaging in a closed cupboard to protect them from dust, moisture, and temperature fluctuations.

CAUTION: Material damage during transport: If not packaged correctly, the instrument may be damaged during transportation. Use the original packaging to avoid damage to the instruments during transportation.

	Temperature	Rel. humidity
Storage	-20°C ... 60°C	10% ... 90%
Application	10°C ... 40°C	30% ... 100%

10. Transport

NOTE: Use the original packaging for shipping in order to avoid damage.

11. Directive compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible office.

12. Disposal

To dispose of these medical devices, no special measures are necessary. Observe the national laws and regulations.

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Symbol description

- WARNING:** Failure to observe may result in injury or even death.
- CAUTION:** Failure to observe may result in damage to or even destruction of the product.
- NOTE:** Special information on the operation of the instrument.
- Read accompanying documents before use
- Manufacturer

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- WARNING:** Risk of infection: KARL STORZ instruments are provided non-sterile and must be cleaned and sterilized/high-level disinfected before initial use and before every subsequent use.
- WARNING:** Risk of infection for patient. Prior to sterilization inspect the instrument for visible soiling (blood, tissue remnants), which indicates that cleaning and disinfection have not been performed or have been performed incorrectly. Repeat cleaning as necessary.
- WARNING:** Risk of infection for patient/staff: Ensure the sterility of the instruments and surroundings.
- WARNING:** Risk of infection or damage to products: Effective reprocessing cannot be guaranteed if non-validated procedures are employed.
- WARNING:** Risk of infection from unsterilized instruments. Only use the sterilization methods recommended in the instruction manual.
- WARNING:** Risk of injury to the patient/user and/or damage to the products: Combinations of medical products are only assured to be safe if they are identified as such in the respective instruction manuals or the intended purpose and interface specifications of the products used in combination permit this.
- WARNING:** Risk of injury to the patient/user. Damaged or incomplete instruments/accessories may cause malfunctions. Always check instruments/accessories before and after use to ensure that nothing is missing/damaged. Replace incomplete/damaged instruments/accessories and do not continue to use them. Before every use check that all of the instruments required for the operation are functioning correctly.

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1. Hold the opening aid ⑪ parallel to the take-apart instrument as shown in Fig. F.
2. Push the opening aid over the instrument until the notch of the opening aid lies directly on the lug ③. To facilitate handling, open the grip of the instrument and press the opening aid gently against the edge of the upper sheath half (see Fig. G).

- The locking mechanism (the lug) is now open (see Fig. H).
- 3. The small ring component can now also be removed (see Fig. I).

5 Assembly

1. Re-attach the disassembled ring component to the body of the instrument. The grip clicks audibly into place in the body of the instrument (see Fig. I).
2. There is a guide in the middle of the slide of the instrument and a corresponding recess in the middle of the body of the instrument. Press the guide into the recess (see arrow in the photograph, see Fig. J).
3. Press the slide of the instrument so that the circular notch is located below the lug (see arrow in the photograph, see Fig. K).
4. Press the circular notch of the slide into the lug. This may require moving the grip of the instrument to and fro a little so that the circular notch of the slide can engage in the lug. Ensure that the notch is fully engaged, as otherwise the safety bar cannot be closed in the next step of assembly.
5. Close the safety bar (see arrow direction in Fig. L).

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6.3 Pre-cleaning at the point of use

- 1 Immediately after a procedure, at the point of use (e.g. operating room), wipe the instruments with a soft, lint-free, disposable cloth, moistened with distilled/demineralized water or cleaning detergent to remove all heavy soiling.
- 2 Contain and transport the contaminated instruments to the decontamination area for manual cleaning.
- 3 If immediate cleaning is not feasible, place the instruments in a container and soak in neutral/ mild pH enzymatic cleaning solution (e.g. Enzol), diluted with distilled/demineralized water to proper concentration per manufacturer's instructions

6.4 Manual cleaning

1. Disassemble the instruments
2. Thoroughly rinse instruments with distilled/ demineralized water to remove gross debris.
3. Completely immerse instruments in a mild/neutral pH enzymatic solution. Keep immersed for 10 minutes.
4. While immersed, brush the instrument with the cleaning brush
5. After immersion period, rinse the instruments with distilled/demineralized water.
6. Dry the instruments with a lint-free soft cloth or compressed air.

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
- ① Jaws
- ② Working insert
- ③ Sheath tube
- ④ LUER-Lock port for irrigation
- ⑤ Cap for irrigation port (Art. No. 8240990 or 29100: pack of 10)
- ⑥ Wheel for turning the jaws
- ⑦ Locking button for releasing the working insert
- ⑧ HF connection
- ⑨ Handle

1 Intended use

CLICKLINE scissors and forceps are used to dissect and grasp or cut, CLICKLINE curettes for dissection and preparation of soft tissue during endoscopic interventions. CLICKLINE scissors and forceps with an HF connection are also used for monopolar coagulation of tissue or for coagulation of minor hemorrhages

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2 Contraindications

Contraindications relating directly to the product are currently not known. Use of CLICKLINE scissors, forceps and curettes is considered to be contraindicated if such application would endanger the patient, e.g. due to the patient's general condition, or if the endoscopic method as such is contraindicated. CLICKLINE scissors and forceps must not be used for interventions on the CNS.

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3 Supplied contents

CLICKLINE instruments each consist of a handle and outer sheath with working insert (instruments in 2, 3 and 3.5 mm sizes, curved instruments) or handle, outer sheath and working insert (5 mm and 10 mm instruments, straight and offset).

4 Application

KARL STORZ recommends the use of the CLICKLINE grasping forceps and scissors with HF instruments from the AUTOCON® range from KARL STORZ.

A peak voltage of 3 KvP must not be exceeded in the Forced Coagulation and Spray Coagulation modes, as this can damage the instrument.

5 Safety information

WARNING: Risk of injury and damage to products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of injury: There is a risk of injury if the insulation of HF instruments is damaged or missing or if there is unintentional contact between the application part and the patient. Do not use HF instruments with damaged or missing insulation for bipolar resection and make sure that the application parts do not come into contact with conductive instruments and accessories. Never place the application parts on the patient.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional raw surfaces, sharp corners, burled edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination indicates that preparation has not been performed, or has been performed incorrectly. Prepare the instruments before initial use and before and after every subsequent use using validated procedures.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

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CLICKLINE 器件

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury and risk of damaging the products: If application parts are used outside of the field of vision there is a risk that tissue and accessories could be damaged unintentionally. Always hold the application parts of the active electrode, laser and other instruments which transmit energy in a target-oriented manner in the field of vision during application.

WARNING: Risk of injury and risk of damaging the products: Sparks may be produced when HF devices and accessories are used. These can cause combustible or flammable liquids to ignite or explode. Make sure that combustible or flammable agents have evaporated before using HF instruments and that there are no flammable gases.

CAUTION: CLICKLINE instruments with HF connection used in combination with HF instruments are intended for a recurring peak voltage of max. 3 kVp and are only suitable for the short-term coagulation of small hemorrhages. Applications with device settings above 3 kVp can damage the instrument. (In acc. with IEC 60601-2-2, 3rd edition). CLICKLINE instruments with HF connection models 50220xx, 50221xx, 50222xx, 50235xx, 50236xx, 50510xx are intended for a recurring peak voltage of max. 2.0 kVp.

CLICKLINE instruments with HF connection models 33125/33125 P, 33127/33127 P are intended for a recurring peak voltage of max. 1.8 kVp.

WARNING: Always use the lowest possible current setting at which adequate cutting and coagulation can be achieved.

WARNING: Tissue desiccation/coagulation can effectively be achieved with either „CUT“ (low voltage) outputs or „COAG“ (high voltage) outputs. „COAG“ outputs are generally associated with much higher open-circuit peak-to-peak voltage than the „CUT“ mode at identical outputs.

WARNING: Use of „cut“ mode only is strongly recommended to reduce capacitive coupling and unintentional tissue effects including those remote from the operative field.

WARNING: If the tissue effect is not being created at nominal settings, do not increase power output until ensuring that the cord and electrodes are in good condition and properly connected and that the dispersive electrode is in proper position and in contact with the patient.

WARNING: It is strongly recommended that the surgeon use high quality surgical gloves that have good insulation properties, thereby minimizing the risk of puncture and burn.

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CLICKLINE 器件

WARNING: Continuous activation of the active electrode while it is not in contact with the tissue may increase the possibility for capacitive coupling, instrument breakdown, and undesirable and remote tissue effects. This is especially true for „COAG“ (high voltage) outputs. Avoid continuous activation of the electrode to reduce the risk of current leakage and diversion to alternate sites.

WARNING: Gas emboli and cardiac arrest have been reported during electrosurgical procedures performed in conjunction with gas insufflation. Careful monitoring of gas volume and patient pressure should be performed throughout electrosurgical procedures in order to minimize any such risk. KARL STORZ recommends that gas insufflation NOT be used during operative hysteroscopic procedures. In addition, minimizing instrument exchanges and limiting the use of the Trendelenberg position may also reduce the risk of gas emboli.

WARNING: Activation of electrosurgical unit simultaneously with the aspiration/irrigation mode may increase the capacitive coupling through the shaft area of the insulated instrument, and cause arcing at the distal end of the instrument.

WARNING: Use of electrosurgical devices can cause severe electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions should be taken to ensure that the patient's well-being is maintained in the event of such interference.

WARNING: If a reduction in the expected electrosurgical output is noted, check the return plate (patient plate or pad, dispersive electrode) to ensure proper positioning. Check the active electrode insulation for possible damage before increasing the power of the electrosurgical generator.

WARNING: The patient should not come into contact with any metal parts that are grounded.

WARNING: Electrodes and probes of monitoring, stimulating, and imaging devices can provide paths for high-frequency currents even if they are battery powered, insulated or isolated at 60 Hertz. The risk of burns can be reduced but not eliminated by placing the electrodes as far away as possible from the dispersive electrode (e.g. patient return electrode, patient plate or pad).

6 Disassembly

1. Open handle completely. If necessary, unlock the locking mechanism on the handle.
2. Press and hold the locking button ⑦ (do not turn).
3. The handle is now opened completely in order to push the sheath tube ③ out towards the front (see Fig. A).
4. To remove the inner section, hold the jaws ① tightly and open the bayonet mount by rotating the sheath tube ③ ¼ turn.
5. Pull out the inner section (see Fig. B)

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
7

Preparation

CAUTION: When preparing and using the solutions, follow the chemical manufacturers' instructions, paying close attention to concentration and exposure times. Prolonged immersion and incorrect concentration may result in damage. Regarding the microbiological effectiveness please ask the manufacturer.

NOTE: The instruments must be completely dismantled prior to preparation.

NOTE: When carrying out any work on contaminated instruments, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: For efficient cleaning, we recommend an ultrasound bath for a minimum of five minutes at 35 kHz and 40 °C.

NOTE: Only process the instruments through manual cleaning OR machine cleaning. Processing through both methods is not necessary.

7.1 Manual cleaning

7.1.1 Pre-cleaning at point of use

1. Immediately after a procedure, at the point of use (e.g. operating room), wipe the instruments with a soft, lint-free, disposable cloth, moistened with distilled/ demineralized water or cleaning detergent to remove all heavy soiling.
2. Contain and transport the contaminated instruments to the decontamination area for manual cleaning.
3. If immediate cleaning is not feasible, place the instruments in a container and soak in neutral/mild pH enzymatic cleaning solution (e.g. Enzol), diluted with distilled/demineralized water to proper concentration per manufacturer's instructions.

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CLICKLINE 器件

7.1.2 Cleaning

1. Disassemble the instruments
2. Thoroughly rinse instruments with distilled/ demineralized water to remove gross debris.
3. Completely immerse instruments with jaws open in a mild/neutral pH enzymatic solution. Flush the lumen with cleaning solution and keep immersed for 10 minutes.
4. While immersed, brush the instruments with a cleaning brush. Brush the lumen with a small brush until no debris can be seen on the bristles of the brush.
5. After immersion period, rinse the instruments with distilled/demineralized water. Flush the lumen with water to remove the cleaning solution.
6. Dry the instruments with a lint-free soft cloth or compressed air. Flush the lumen with compressed air to dry them out.

8 Sterilization

WARNING: KARL STORZ recommends that the instruments be sterilized disassembled. However, Clickline instruments may be sterilized assembled in a steam pre-vacuum sterilizer only.

The following steam sterilization parameters have been validated by KARL STORZ for use in the U.S.:

Pre-Vacuum

- Temperature: 132 °C
- Exposure Time: 4 minutes

Gravity Displacement

- Temperature: 121 °C
- Exposure Time: 45 minutes

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


CLICKLINE 器件

9 Assembly

NOTE: Safety feature! For safety reasons, handles with an HF connection can only be connected to insulated outer sheaths.

1. Insert the inner section (2) from the front into the sheath tube (3) (see Fig. D).
2. Lock the inner section (2) by rotating ¼ turn in the bayonet mount on the sheath tube (3) (see Fig. E).
3. Insert the sheath with the jaws closed as far as it will go into the completely opened handle (see Fig. F).
4. Then lock the sheath tube by closing the handle completely (see Fig. F).
5. The handle must be lubricated regularly with care oil (# 27656 B), in order to ensure that the handle moves smoothly (see Fig. G and H).
6. Also apply a drop of care oil (# 27656 B) to the working insert of the instrument.
7. Test for proper operation.

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

CLICKLINE 器件

10 Accessories

The springs 33131 S ⑩ and 33131 F ⑪ can be attached to the handle 33131.

CAUTION: The long leg ⑪ of the springs 33131 S and 33131 F must point toward the distal end of the assembled CLICKLINE instrument. Attached thus, the spring does not get caught in the fulcrum of the handle arms when closing the handle.

33131 S is for opening the jaws, while 33131 F automatically closes the jaws by spring tension (see Fig. J and K). The wheel 33120 ⑬ (pack of 5) can be fitted on any CLICKLINE handle. The wheel makes rotating the outer sheath easier (see Fig. L).

Recommended unipolar high frequency cords	
26002 M	Unipolar High Frequency Cord, with 4 mm plug, length 300 cm, for models KARL STORZ, Erbe type T, older models and Ellman
26004 M	Unipolar High Frequency Cord, with 4 mm plug, length 300 cm, for use with Martin HF units
26005 M	Unipolar High Frequency cable, with 5 mm plug, length 300 cm, for AUTOCON® II 400 SCB system (111, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ AUTOCON® system (50, 200, 350) and Erbe type ICC
26006 M	Unipolar High Frequency Cord, with 8 mm plug, length 300 cm, for use with HF units, models KARL STORZ AUTOCON® II 400 SCB (112, 116) and Valleylab

11 Application/Storage conditions

	Temperature	Rel. humidity
Storage	-20 °C ... +60 °C	10 % ... 90 %
Application	+10 °C ... +40 °C	30 % ... 70 %
Air pressure	700 hPa – 1060 hPa	

12 Disposal

To dispose of CLICKLINE instruments, no special measures are necessary. National laws and regulations must be observed.

13 Directive compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

14 Packaging symbols

For the meanings of the symbols printed on the packaging, please refer to the "Packaging symbols" accompanying instruction leaflet 96216316DF. This can be downloaded from www.karlstorz.com.

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Normal use

The EndoCAMEleon® is a special **HOPKINS®** telescope which, thanks to a swiveling V-block (3), allows the field of view to be set between 0° and 120° during endoscopic examinations and interventions. The EndoCAMEleon® is used, amongst other things, for illuminating and observing the operating field in laparoscopy, gynecology, urology, pediatrics and thoracoscopy.

NOTE: It is recommended to use the EndoCAMEleon® in combination with the fiberoptic light cable 495NCS (4.8 mm).

No **contraindications** relating directly to the product are currently known. Use of the instrument is contraindicated if, in the opinion of the responsible physician, such an application would cause danger to the patient, e.g., due to the general condition of the patient, or if the method as such is contraindicated. Unauthorized conversions or modifications to the instrument are prohibited for safety reasons.

Explanation of the arrow icons (4)

The following fields of view result when you rotate the setting ring (2) to one of the following arrows:

- ↑ Corresponds to the field of view of a 0° telescope
- ↖ Corresponds to the field of view of a 30° telescope
- ↗ Corresponds to the field of view of a 45° telescope
- ↘ Corresponds to the field of view of a 70° telescope
- Corresponds to the field of view of a 90° telescope
- ↙ Corresponds to the field of view of a 120° telescope

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Inspection

Inspect the telescope for damage immediately after receipt and before and after each use.

Image waveguides

Examine the optical end faces (distal [=closest to the patient] tip and eyepiece) for scratches and surgical or cleaning residues. The faces should be smooth and shiny.

In order to examine the image quality, the telescope should be rotated slowly whilst being looked through. If the image is completely or partially impaired, either the rod lens system in the jacket tube is broken or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

Light cables

Hold the light cable connection **by the proximal (=furthest from the patient) end** in the direction of a light source and inspect the number of black dots on the light outlet points of **the distal (=closest to the patient) tip**. Such dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image.

WARNING: If half of the optical fibers of a light outlet point are broken then this results in a significant deterioration in the brightness and thus the image quality in this area. For example, the surgeon should not operate using an EndoCAMEleon® in which the rearmost light outlet point is completely broken.

WARNING: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy or if no image or only parts thereof are visible.

Sheath

Inspect the entire sheath for dents and scratches. The sheath must not be bent.

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Application

WARNING: The instruments may only be used by persons with an appropriate medical qualification and who are acquainted with the application.

WARNING: The product described here must be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING: The instruction manuals and interface specifications for the medical devices used in combination must be observed precisely.

WARNING: Infection risk for patients/personnel: These instruments are **not sterile** when delivered. Instruments must be cleaned, disinfected and sterilized by validated procedures before they are used for the first time, and before and immediately after every subsequent use.

WARNING: Risk of infection to patients: Check the instrument for any visible soiling (blood, tissue residue), which indicates preparation was not carried out or was carried out incorrectly. If necessary, carry out preparation again (cleaning, disinfection, sterilization).

WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (=furthest from the patient) of the telescope is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the setting ring ② which can cause cardiac pacemakers or other electronic devices to malfunction.

WARNING: Risk of injury for patients/users: Never use excess force on the instrument. Bent instruments must be replaced and not bent back into shape.

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WARNING: Risk of injury for patients/users: Do not use damaged instruments/accessories (e.g., if they have rough surfaces, sharp corners, burred edges, or projecting parts). Always check instruments/accessories before and after use to ensure that nothing is missing. Test for proper operation. Replace incomplete instruments/accessories and do not reuse them.

WARNING: Risk of injury for patients: Prior to every use, check that the components are securely assembled by gently pulling and/or turning in both directions without applying too much force.

NOTE: Every time the EndoCAMEleon® is used, it should first be newly adjusted by turning the setting ring ② once fully clockwise and once fully counter-clockwise.

NOTE: Check the functioning of the EndoCAMEleon® prior to use by turning the setting ring ②. It must be possible to steplessly adjust the field of view between 0° and 120°, where applicable by clicking to the different directions of view with minimal force.

WARNING: Risk of injury for patients: The high level of light energy means that the distal end, the light port and adjacent components may become so hot that they can cause skin burns to the patient or user! Avoid tissue coming into direct contact with the distal end and/or light port of the telescope.

WARNING: Risk of injury for patients: Thermal tissue damage: The higher the intensity of the light source, the greater the heat emission. For this reason, the light source output should always be selected to be as low and short as is possible while still ensuring optimal illumination of the operating field.

WARNING: Danger of glare! Risk of burns! Never look into the open end of a connected light cable or endoscope. Both light and heat are also radiated in directions which are not directly in the user's field of vision.

WARNING: It is important that you are always aware of the currently set field of view in order to avoid unintentional tissue damage.

CAUTION: The application and storage conditions must be observed.

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Precautions

CAUTION: In order to avoid breaks, telescopes must always be held by the eyepiece and never by the distal (=closest to the patient) end alone.

CAUTION: Never bend the jacket tube. This can result in breaks or cracks in the rod lens system.

CAUTION: Place the telescopes down with care. Hard impacts, particularly to the distal (=closest to the patient) end, but also to the proximal (=furthest from the patient) end, can cause damage or cracks. This results in cloudy and/or blurred image areas.

CAUTION: To prevent damage to the telescope the setting ring must be turned back to the 0° direction of view after every use. This applies in particular to the preparation and transport of the telescope

CAUTION: Therefore, during preparation, cleaning, disinfection and sterilization, it must be ensured that the telescopes cannot be damaged by other instruments. HOPKINS® telescopes are best placed individually or use containers in which the telescopes can be secured.

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CAUTION: The manual 'Cleaning, sterilization and care of KARL STORZ instruments' must be followed. The cleaning, disinfection and sterilization procedures are explained therein in detail. The up-to-date list of KARL STORZ approved materials can be found on the internet at 'www.karlstorz.com' or in the manual 'Cleaning and sterilization of KARL STORZ instruments'.

Please also observe the ready-reference guide 'Cleaning, disinfection and sterilization of HOPKINS® telescopes/miniature endoscopes'.

CAUTION: When preparing and using the solutions, follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged exposure may result in corrosion and/or material damage.

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Manual cleaning

CAUTION: Telescopes must not be cleaned in an ultrasound bath.

CAUTION: Telescopes must not be placed in physiological saline solution under any circumstances.

Carefully clean the surfaces with a disinfected sponge or a soft brush.

Only cleaned and disinfected instruments can be sterilized effectively.

Immerse all the parts in a disinfectant solution making sure that all the air bubbles have been removed and that the parts are completely covered and taking the chemical manufacturers' instructions on concentration and exposure times into account. Then rinse out thoroughly with completely demineralized water and dry with sterile compressed air.

Machine cleaning

Only use special methods which have been tested and approved (maximum disinfection temperature 93°C). Preferably use a cleaner with a pH in the neutral range (e.g., an enzymatic cleaner). Ensure that all parts are stored securely to avoid any damage occurring when machine preparation is in progress. To this end, the parts should ideally be placed in a tray holder.

The use of the following processes is permitted:

- chemical disinfection
- thermal cleaning.

Sterilization

These instruments can be steam sterilized at up to 134°C ± 3°C. The instrument can be sterilized in its assembled state using the fractionated vacuum method. The instructions for use provided by the sterilization device and disinfectant manufacturers must be followed exactly.

When the sterilization cycle is complete, allow the container to cool down to room temperature.

CAUTION: If the cooling phase is accelerated (e.g., using compressed air), the telescope may be damaged. Never pour cold water over the telescope.

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NOTE: If deactivation of Creutzfeldt-Jakob-disease (CJD) is necessary, the sterilization method recommended by WHO (WHO/CDS/PH/92.104) may be used (18 minutes at 134 °C). Extension of the sterilization time up to one hour at 134 °C is possible in special cases.

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Application/Storage conditions

	Temperature	Rel. humidity
Storage	-20°C to +60°C	10% to 90%
Application	+10°C to 40°C	30% to 70%

Directive compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

Disposal

To dispose of this medical device, no special measures are necessary. National laws and regulations must be observed.

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1 Intended use

The fixation button is used for reconstructive treatment of ligament ruptures as well as for femoral fixation of an implanted cruciate ligament replacement.

Indications

Reconstructive therapy of ruptures to the anterior and posterior cruciate ligament by means of autologous grafts as well as reconstructive treatment of the AC joint.

Contraindications

- Osteoporosis and osteomalacia
- Degenerative osteopathies
- Adiposity or patient overweight, leading to excessive strain on the fixation button
- Osteomata in the area in which the fixation button is to be placed
- Deformities of the bone, or general conditions of the bones which exclude implantation of a fixation button in the opinion of a physician
- Systemic diseases and metabolic disorders

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2 Factors that may compromise the success of the operation

- Inadequate bone quality, (e.g. osteoporosis, previous cruciate ligament operation)
- Systemic diseases or metabolic disorders
- History of infections or falls
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Local bone tumors

Possible adverse effects

- I. Infection, both deep and superficial
- II. Allergies, mild inflammatory and foreign body reactions to implant material

Impairments to imaging examinations

The quality of the image produced when conducting certain imaging examinations, such as MRI or CT, can be impaired by the implants.

The patient must be made aware that, before undergoing such an examination, he should inform the doctor responsible that he has an implant of this type.

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3 Patient information

The physician must inform the patient of the factors which may compromise the success of the operation and possible adverse effects as listed under items "Contraindications" and "Possible adverse effects" as far as they apply to the individual patient. The patient should also be informed about any measures he can take to minimize said factors and adverse effects and the possible consequences.

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4 Surgical technique

Surgeons who are not familiar with the specific product and instrument set are expected to study the surgical technique thoroughly before performing the operation. The product-specific documentation is available from the responsible KARL STORZ representative.

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5 Handling

 **NOTE:** This product is packaged in protective packaging and is sterile.

 **WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

 **WARNING:** Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

 **WARNING:** Before use, check the expiration date and ensure that the packaging is undamaged. Products which have expired or which are damaged must not be used.

 **WARNING:** Risk of infection: The use of non-sterile instruments poses a risk of infection for patients. Do not remove the products from the sterile packaging until immediately before use.

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 **WARNING:** Risk of infection. Reprocessing of disposable articles can jeopardize the functionality and safety of the product and lead to infections. Disposable articles must not be reprocessed. Dispose of disposable articles properly.

 **WARNING:** Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional raw surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

 **WARNING:** Risk of injury: Overloading the instrument by exerting too much force on it may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

 **NOTE:** The tibial fixation button is only used to assist the normal healing process. It is not capable of bearing the complete load of bodily weight before the cruciate ligament replacement graft has completely settled. The patient must be instructed accordingly.

6 Application

During implantation of the extracortical fixation button, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on progress, as well as monitoring the vital signs and the anesthetic.

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7 Traceability

The self-adhesive identification labels included in the original packaging must be retained for purposes of traceability; they should be placed in the operation report and the patient's medical record.

The individual responsible for the stockroom in the hospital must keep a record of the available products, incl. article number (REF) and batch number (LOT) to ensure that fast action can be taken in the event of changes or product recalls.

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8 Product Material

The extracortical fixation button is made of a high-quality, implantable titanium alloy (Ti 6 AL 4V or F 136 ASTM)

 **NOTE:** It is a sterile disposable product.

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9 Packaging and sterility

Each individual product is double packed and sterilized using either gas (EO) or gamma rays with a dose of at least 25 kGy. To find out which procedure has been used in each case, please refer to the symbol on the packaging. Sterility is guaranteed up to the expiry date, provided that the packaging has not been opened or damaged.

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11 Storage conditions

The product must be stored in the original unopened packaging in a dry place at room temperature, and after the expiration date it must be resterilized before use.

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12 Liability

Complications or other adverse affects which might occur due to incorrect indication or unsatisfactory surgical technique, selection of unsuitable products or instruments or handling of them, asepsis, etc. come under the liability of the surgeon.

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Karl Storz GmbH & Co. KG