Notariat Tuttlingen B2

Amtsverwalter Dennis Müller

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1 Important information for users of KARL STORZ devices	1
1.1 Revision index	1
1.2 Validity	1
1.3 Other applicable documents	1
1.4 Symbols and labeling	2
1.4.1 Structure of warnings	2
1.4.2 Risk levels in warnings	2
1.4.3 Tips	2
1.4.4 Other symbols and labeling	3
2 Safety	4
2.1 Intended use	4
2.2 Contraindications	4
2.3 General safety information	5
2.4 Personal safety instructions	7
2.4.1 Ambient conditions	7
2.4.2 Patients with pacemakers	7
2.4.3 Hazard-free patient positioning	7
2.4.4 Correct connection of the HF device	8
2.4.5 Correct use of the HF device	9
2.4.6 Configuring HF device settings and using accessories	9
2.5 Product-related safety instructions	10
2.6 Safe handling (general instructions)	11
2.6.1 Operation area: Avoiding ignition and explosions	12
2.6.2 Applying the neutral electrode	12
3 Description	16
3.1 Display and control elements	16
3.1.1 Control elements on the front panel	16
3.1.2 Unipolar socket module (left)	17
3.1.3 Bipolar socket module (right)	17
3.1.4 Control elements on rear panel	18
3.2 Symbols used on the device	19
3.2.1 Type plate	21
3.2.2 Symbols on the packaging	21
3.3 Scope of supply	21
3.4 Components required for operation	21
3.5 Operating conditions	21

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4 Preparation	22
4.1 Setting up the HF device	22
4.2 Switching on the HF device	23
4.3 Connecting instruments	24
4.3.1 Instruments for unipolar application	24
4.3.2 Instruments for bipolar application	25
4.3.3 Connecting the footswitch	25
4.4 Test for proper functioning	26
4.4.1 Automatic test function	26
4.4.2 Performing the test for proper functioning	26
4.4.3 Actions in case of problems	27
4.5 Neutral electrode monitoring	27
4.5.1 General information	27
4.5.2 EASY neutral electrode monitoring (EASY monitoring)	28
5 Operation	29
5.1 Connecting the power supply	29
5.2 Program overview	29
5.2.1 Display	29
5.2.2 Status bar	30
5.3 Activating and deactivating sockets	30
5.4 Unlocking the screen	31
5.5 Configuring output currents	32
5.5.1 Selecting the mode	32
5.5.2 Specifying power limits	33
5.5.3 Selecting the effect	33
5.5.4 Assigning the footswitch	34
5.5.5 Selecting the neutral electrode	36
5.5.6 Coding system	38
5.6 Mode overview	39
5.6.1 Unipolar modes	39
5.6.2 Bipolar modes	42
5.7 Unipolar cutting modes	44
5.7.1 Laparoscopy	44
5.7.2 Standard	44
5.7.3 Micro	45
5.7.4 Resection	45
5.7.5 Resection C-Cut®	46
5.7.6 SupraLOOP	46
5.7.7 Dry	47
5.7.8 Gastro Loop 1	47
5.7.9 Gastro Loop 2	48

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Important information for users of KARL STORZ devices

1 Important information for users of KARL STORZ devices

It is recommended to check the suitability of the products for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual forms part of the device. KARL STORZ GmbH & Co. KG, referred to in the following simply as KARL STORZ, assumes no liability nor provides any warranty whatsoever for damage and consequential damages that arise due to non-compliance with the instruction manual.

- Read the instruction manual carefully before using this device.
- Store the instruction manual in a safe place throughout the service life of the device.
- Keep the instruction manual accessible to surgical staff.
- Give the instruction manual to each successive owner and/or user of this device.

1.1 Revision index

Software version	Last revised
Valid from version 2.2.0	Version 1.0 – 2016/04

1.2 Validity

This instruction manual Art. no. 96206584 DER applies only to the devices UH 400, UH 400U, UH 401 and UH 401U.

1.3 Other applicable documents

- Comply with other applicable documents in the appendix.

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1.4 Symbols and labeling

1.4.1 Structure of warnings SIGNAL WORD

Risk type, source and consequences thereof (Personal injury)!

- Measure to avoid the risk

NOTE

Risk type, source and consequences thereof (Material damage)!

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Symbol	Risk level	Probability of occurrence	Consequences of non-compliance
	DANGER	Immediate risk	Death, serious injuries
	WARNING	Possible risk	Death, serious injuries
	CAUTION	Possible risk	Minor injuries
	NOTE	Possible risk	Material damage

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1.4.3 Tips

- Tips and additional information to facilitate tasks

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**Important information
for users of KARL STORZ
devices**

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1.4.4 Other symbols and labeling

Symbol/ labeling	Meaning
☑	Prerequisite for an action
▷	Action with one step
1. 2. 3.	Action with several steps in a binding sequence
☛	Result of preceding action
•	List (first level)
•	List (second level)
Emphasis	Emphasis
..., see section xxx page xxx	Cross reference

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2 Safety

2.1 Intended use

The HF device is intended exclusively for the generation of electrical power for unipolar and bipolar cutting and coagulation on tissue structures in surgical operations.

It is used in the following areas:

- General surgery
- Arthroscopy
- Endoscopy
- Gastroenterology
- Gynecology
- Hand surgery
- ENT
- Cardiac surgery (including open-heart surgery)
- Minimally invasive surgery (MIS)
- Neurosurgery
- Orthopedics
- Plastic surgery and dermatology
- Pneumology
- Thoracic surgery
- Urology, including transurethral resection (TUR)
- Oral and maxillofacial surgery

2.2 Contraindications

Not intended for the preparation of metallic implants such as stents.

- Not suitable for surgical procedures which require continuous activation of the HF generator.

The following body parts and tissue types must not be treated:

- Not suitable for working on bone tissue.

The HF device must not be used if, in the opinion of an experienced physician or according to current specialist literature, such an application would endanger the patient, e.g. due to the patient's general condition, or if other contraindications are present.

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-  KARL STORZ requires that the HF device is operated under the supervision of qualified and authorized personnel. The surgeon and medical staff must be trained in the fundamental principles, rules for use and risks of HF surgery and must be familiar with these in order to safely and reliably prevent putting patients, staff and equipment at risk.
-  Contact your KARL STORZ specialist dealer about training and training material.
-  Any other use is neither intended nor proper and must be effectively prevented.

2.3 General safety information

- ▷ Ensure that no electronic devices that are subject to interference from electromagnetic fields are set up in the vicinity of the HF device.
- ▷ Observe the instructions on electromagnetic compatibility provided in section EMC, page 221.

WARNING: Danger of electric shock!

- ▷ Always connect the HF device to a mains power system with a protective conductor in order to prevent electric shock.

Additional equipment connected to medical electrical equipment must verifiably comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the

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
above-mentioned standard requirements. If in doubt, please contact your local specialist dealer or the Technical Support, see section Technical Support, page 100.

-  To protect personnel, KARL STORZ recommends the use of a smoke evacuator to extract electrosurgical smoke. To this end we recommend the KARL STORZ S-PILOT™ UP 501 which, in combination with the AUTOCON® III 400 (plus connecting cable UP 004), guarantees automatic smoke extraction.

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2.4 Personal safety instructions

2.4.1 Ambient conditions

- ▷ Do not use the HF device in the immediate vicinity of the patient. Observe the minimum distances recommended by KARL STORZ (see Fig. 1).

2.4.2 Patients with pacemakers

Malfunction or destruction of the pacemaker can endanger the life of the patient or result in irreversible injuries to the patient.

- ▷ In the case of patients with pacemakers, consult the cardiologist before carrying out HF surgery.
- ▷ Use bipolar HF methods.
- ▷ Attach the HF neutral electrodes close to the operating field.
- ▷ Set the demand pacemaker to a fixed frequency.
- ▷ Ensure that the pacemaker does not come into contact with the HF electrode.
- ▷ Keep a fully operational defibrillator within reach.
- ▷ Carry out a postoperative pacemaker check.

2.4.3 Hazard-free patient positioning

- ▷ Position patients so that they are not touching any metal parts that are grounded or have considerable capacitance relative to ground (e.g. operating table brackets). If necessary, place anti-static cloths between the patient and the underlay.
- ▷ Ensure that the patient does not touch any wet clothes or underlay.
- ▷ Place anti-static cloths between areas with heavy sweating and skin-to-skin contact areas on the patient's torso.
- ▷ Ensure that the patient is resting on a suitable surface in order to prevent pressure necrosis.
- ▷ Drain urine via a catheter.

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2.4.4 Correct connection of the HF device

- ▷ Always ground the HF device via the potential equalization connector. Also observe the requirements in section 8.6.7 of IEC 60601-1 regarding medical electrical systems.

i The device features a potential equalization connector which must be connected via a potential equalization line with a connector in the room where it is set up. This additional potential equalization measure means that all conductive parts in the vicinity of the patient are connected electrically to each other and small contact voltages are compensated. As such, both transient disturbances from the switching of adjacent devices and the risk to patients from low-frequency electric current are reduced. The connection should be as low-resistance as possible and should not exceed a resistance value of 0.2 Ohm between the central potential equalization of the operating room and the potential equalization connector AUTOCON® III 400 and every other metallic touchable surface. If the design of the installation does not include an equipotential bonding bar, the connection to the PE conductor should be as low-resistance as possible.

- ▷ With simultaneous use of HF surgical devices and physiological monitoring devices on a patient, the monitoring electrodes should be positioned as far away as possible from the surgical electrodes.
- ▷ Needle electrodes for monitoring are not recommended.
- ▷ In all cases monitoring systems are recommended which include mechanisms for limiting the high-frequency current such as protective resistors and HF chokes.
- ▷ Place lines from monitoring devices so that they do not lie on the patient's skin.
- ▷ Keep the leads to the HF electrodes as short as possible and position them so that they do not touch the patient or other leads.
- ▷ Do not place any objects on the HF device.

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2.4.5 Correct use of the HF device

Inadvertent activation of the HF device outside the user's field of vision can injure the patient.

- ▷ Activate the HF device only when the electrode is in your field of vision and you can quickly deactivate the HF device at all times.
- ▷ If the HF device is activated inadvertently, switch it off immediately using the Standby button.
- ▷ Take particular care when using a footswitch or manual switch.

Improper preparation or errors in the HF device can cause damage to the HF device.

- ▷ Use the automatic monitoring functions to ensure that the HF device is working properly. See the section Test for proper functioning, page 26 for information on the automatic test functions.
- ▷ Ensure that no conductive fluids (e.g. blood or amniotic fluid) have penetrated the footswitch or the manual switch.
- ▷ Ensure that the cables for the footswitch and the manual switch are free from short circuits and broken leads.

2.4.6 Configuring HF device settings and using accessories

Setting the output power too high can injure the patient! Therefore, before you increase the output power, ensure that:

- the neutral electrode is attached properly;
- the working electrodes are clean;
- the plug connections are all correct.

A clearly lower output value or functional failure of the HF surgical unit during normal operation can be due to insufficient attachment of the neutral electrode or insufficient contact in the connections. In this case, check attachment of the neutral electrode and the connections before selecting a higher output power.

Setting the HF device correctly

- ▷ To prevent inadvertent (thermal) tissue damage during operations on body parts with small cross sections and in areas with high resistance (bones or joints), use the bipolar method in these areas.

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- ▷ Set the level of the acoustic signal that sounds when the electrode is activated so that it is always clearly audible.

Nerve and muscle excitation due to low-frequency currents!

During HF surgical operations (especially when an arc is formed), part of the HF current is converted into a low-frequency current. This current can trigger muscle contractions in the patient:

- ▷ To minimize the risk of injury to the patient, set the power and the effect as low as possible.

Correct use of accessories

- ▷ Only use insulated accessories.
- ▷ Check all electrodes for sharp edges and projecting parts before use.
- ▷ Only use electrodes that are free of defects and in good working order.
- ▷ Never place active electrodes on or near the patient.
- ▷ Do not remove hot electrodes from the patient's body directly after cutting or coagulation.
- ▷ Ensure that there is sufficient distance between the patient cables and the cables of the HF device.
- ▷ Do not run the patient cables across the patient.

2.5 Product-related safety instructions

Devices manufactured by KARL STORZ are developed in accordance with the current state of technology and generally accepted safety rules. Despite this, using these products can lead to risks to the life and health of the user or third parties and/or damage to the device or other objects.

- ▷ Only use accessories approved by KARL STORZ, see section Accessories/Replacement parts, page 216.
- ▷ Only use the device if it is free of technical defects and in good working order and only for the intended purpose, always remaining aware of safety requirements and risks and complying with this instruction manual.

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▷ Have faults that can adversely affect safety (e.g. deviations from the permissible operating conditions) rectified without delay.

▷ Wipe down the HF device only with cleaning agents and disinfectants that are approved in the country of use for surface cleaning. Other reprocessing methods are not permitted! (See section Cleaning and disinfection, page 96.)

▷ If any fluids penetrate the device, drain them immediately.

Damage to the device can lead to an undesirable increase in output power due to a device malfunction.

Certain devices or accessories can cause danger in lower power settings. For example, the risk of gas embolism in argon assisted coagulation rises if the HF power is insufficient for the fast creation of an impenetrable eschar layer on the target tissue.

2.6 Safe handling (general instructions)

- ▷ Before each use of the device, check to ensure that it is functioning properly and is in good working order and connected properly.
- ▷ Observe the instructions on intended use in conformance with standards, see section Error indications for EASY monitoring, page 93.
- ▷ During use, always observe and comply with the acoustic signals and/or error indications of the HF device, see section Error indications for EASY monitoring, page 93.
- ▷ The device and accessories may be operated and used only by people who have the necessary training, knowledge and experience.
- ▷ Regularly inspect the accessories, especially the electrode cables, endoscopic accessories and neutral electrodes, for damage to the insulation, proper operation and expiration date.
- ▷ Never place the instruments on the patient or on the devices.

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- ▷ Ensure that no instruments are being cleaned when AUTOSTART is activated.
- ▷ Wear suitable gloves during operations.

2.6.1 Operation area: Avoiding ignition and explosions

Sparks are generated when the HF device is used as intended!

- ▷ Do not use the HF device in areas where there is a risk of explosion.
- ▷ Do not use any flammable or explosive liquids.
- ▷ If the display fails, do not use the HF device any longer!
- ▷ During operations (e.g. in the head or thoracic regions), avoid using ignitable anesthetics and gases that support combustion (e.g. nitrous oxide or oxygen) or extract them using a vacuum system.
- ▷ Use exclusively non-flammable cleaning agents, disinfectants and solvents (for adhesives). If you cannot avoid the use of flammable cleaning agents, disinfectants or solvents, ensure that they have fully evaporated before using HF surgical equipment.
- ▷ Ensure that no flammable liquids collect beneath the patient or in body cavities (e.g. the vagina). Suction or flush body cavities before activating the device.
- ▷ Wipe off all liquids before using the HF device.
- ▷ Ensure that no ignitable endogenous gases are present.
- ▷ Ensure that all materials saturated with oxygen (e.g. cotton or gauze) are kept far enough away from the HF environment that they cannot ignite.

2.6.2 Applying the neutral electrode

- i** Observe the instructions on the use of the neutral electrode in the instruction manual and the information on the package of the neutral electrode.

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Application point for neutral electrode

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In the unipolar HF method, the neutral electrode feeds the current introduced into the patient's body at the surgical site back to the HF device.

▷ To prevent a rise in temperature at the current exit point, the following conditions must be ensured:

- Sufficiently large contact surface between the neutral electrode and the patient's body;
- High electrical conductivity between the neutral electrode and the patient's body.

▷ To prevent the patient being burned by the neutral electrode, you must comply with the following conditions:

- Select the application point for the neutral electrode so that the current paths between the active and neutral electrodes are as short as possible and run longitudinally or diagonally to the patient's body (because muscles are more conductive in the direction of the fibrils, see Fig. 2).
- For surgery in the thoracic region, do not run the current path transversely across the patient's body and ensure that the patient's heart is never in the current path.
- Depending on the surgical site, apply the neutral electrode to the nearest upper arm or thigh if possible, but never closer than 20 cm.
- In the case of self-adhesive disposable electrodes, comply with any further manufacturer instructions regarding the point of application.
- Ensure that the application area is free of scar tissue, bony protuberances, surface hair and ECG electrodes.
- Ensure that there are no implants (e.g. bone nails, bone plates or endoprotheses) in the current path.
- Ensure that no short circuits can occur at the neutral electrode connection.

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- Avoid areas where fluids may collect.
- As far as possible, use split neutral electrodes with a sufficiently large surface area (patient age and max. output power during operation have to be considered).

Before applying the neutral electrode

- ▷ Remove excess body hair.
- ▷ Clean the application site, but do not use any alcohol, since it dries out the skin and increases the transition resistance.
- ▷ If the patient has poor circulation, massage or brush the application site.
- ▷ Attach the neutral electrode over the entire contact surface evenly. Secure reusable neutral electrodes with rubber bands or elastic straps so that they do not loosen when the patient moves. Ensure that the patient's circulation is not impaired (risk of necrosis).
- ▷ Never use wet cloths or conductive pastes.
- ▷ Ensure that no liquids (e.g. cleaning fluids, disinfectants, blood or urine) penetrate between the patient and the neutral electrode.
- ▷ Do not place the neutral electrode under the patient's buttocks or back.
- ▷ Ensure that there are no ECG electrodes in the current path of the HF device.
- ▷ Check the neutral electrode for damage and to ensure that it is working properly.
- ▷ Replace defective accessories immediately.

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Example application using a disposable electrode

- ▷ Remove the protective film and attach the self-adhesive disposable electrode to the patient. Ensure that the long edge of the disposable electrode faces the operation site and the electrode is fully in contact with the skin. This avoids excessive current concentration on the short edge.
- ▷ Using both hands, press the self-adhesive disposable electrode firmly against the patient's skin.
- ▷ Clamp the electrode tab to the electrode cable.
- ▷ After the operation, remove the disposable electrode carefully to avoid skin damage.

With a one-piece neutral electrode

- ▷ Check the one-piece neutral electrode during surgery.
- ▷ Ensure that the one-piece electrode is not blocked at the device.

With a split neutral electrode

- ▷ Apply the split neutral electrode correctly and without any additional objects, as the HF device does not recognize the bridging of the split surfaces with other objects.
- ▷ Ensure that the HF current flows equally to both parts of the split neutral electrode.

 For monitoring the neutral electrode, see section EASY neutral electrode monitoring (EASY monitoring), page 28.

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3 Description


3.1 Display and control elements

3.1.1 Control elements on the front panel

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① While activating an instrument, the activation bar (7 – 10) of the corresponding socket illuminates yellow or blue.

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- ① Standby button (ON – surrounded by a white light)
- ② Symbol 'Standby button'
- ③ Neutral electrode isolated from ground for HF
- ④ Symbol 'CF type applied part with defibrillation protection'
- ⑤ Symbol 'Observe instruction manual'
- ⑥ Touch screen with mode selection buttons
- ⑦ Activation bar upper unipolar socket
- ⑧ Activation bar lower unipolar socket
- ⑨ Activation bar upper bipolar socket
- ⑩ Activation bar lower bipolar socket

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3.1.2 Unipolar socket module (left)

- ⑪ Connection socket for unipolar instruments with hand or footswitch* (activated state – white surrounding light)
- ⑫ Connection socket for unipolar instruments with hand or footswitch* (activated state – white surrounding light)
- ⑬ Connection socket for neutral electrode* (activated state – white surrounding light)

*Applied part type CF according to IEC 60601-1, with defibrillation protection

Unipolar connection socket ⑪ (see Fig. 3)

- Connection socket ⑪ 'top':
- ⑭ 3-pin US type
- ⑮ 5 mm KARL STORZ/Erbe VIO
- ⑯ 4 mm socket (foot-switched)

Unipolar connection socket ⑫ (see Fig. 4)

- Connection socket ⑫ 'bottom':
- ⑭ 3-pin US type
- ⑮ BOVIE (foot-switched)
- ⑯ 4 mm socket (foot-switched)

Connection socket for neutral electrode (see Fig. 5)

3.1.3 Bipolar socket module (right)

- ⑰ Connection socket for bipolar instruments with footswitch, fingerswitch or AUTOSTART* (activated state – white surrounding light)
- ⑱ Connection socket for bipolar instruments with footswitch, fingerswitch or AUTOSTART* (activated state – white surrounding light)

Bipolar connection sockets

- Upper bipolar socket ⑰ (see Fig. 6)
- ⑲ 2-pin US type (28,58 mm)
- ⑳ 1x KARL STORZ/Erbe VIO
- Lower bipolar socket ⑱ (see Fig. 7)
- ⑲ 2-pin US type (28,58 mm)
- ⑳ 2x KARL STORZ/Erbe VIO

* Applied part type CF according to IEC 60601-1, with defibrillation protection

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
① Disconnect patient-related sockets (⑪, ⑫, ⑬, ⑰, and ⑱) from the patient before commencing service and maintenance work, see section 7 Reprocessing and section 8 Maintenance and repair.

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- ① Connection socket 1 for footswitch
- ② Connection socket 2 for footswitch
- ③ Connection for potential equalization
- ④ Power supply for IEC connector
- ⑤ Fiber-optic signal input socket
- ⑥ Fiber-optic signal output socket

Use the following connections only for service and training purposes:

- ⑦ Ethernet connector (only for connection to KARL STORZ OR1™ control NEO)
- ⑧ USB connector
- ⑨ Connector S-PILOT™ UP 501 with UP 004 S-PILOT™ connecting cable, dia. 3.5 mm, 300 cm (max. 25 V/1 A)
- ⑩ UART communication interface
- ⑪ Power switch
- ⑫ Line fuses (see page 102)

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-  The Ethernet interface is only suitable for connection of the OR1™ control NEO for a dedicated network.
-  Only use the USB interface as described in this manual. To this end, only use the KARL STORZ USB-Stick supplied.

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Symbol	Designation
	Footswitch connection
	Neutral electrode isolated from ground for HF
	CF type applied part with defibrillation protection (cardiac floating)
	Alternating current
	Standby
	During activation (of the HF device) HF energy in the radio frequency range 9 kHz to 400 GHz is applied, which produces electromagnetic radiation. The device uses close-range radio frequency identification technology at the HF connection lines (RFID).
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE), see Disposal.
	Labeling of the (active) HF output; Caution: Hazardous voltage.

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
Symbol	Designation
	Manufacturer
	Date of manufacture
	Observe instruction manual
	Potential equalization connector
	Fiber-optic signal input
	Fiber-optic signal output
	Ethernet connector (connection to KARL STORZ OR1™ control NEO system).
	USB connector*
	S-PILOT™ connector* (surgical smoke extraction)
	UART communication interface for service purposes*
	Reference number
	Serial number
	The device satisfies the Medical Device Directive (MDD) 93/42/EEC and 2007/47/EC (NB 0123) as well as the R & TTE Directive 1999/5/EC.
	Observe the instruction manual
	General warning of physiological effect: The device generates physiological effects which may endanger the patient and/or user.
	'Rx only'

* In the event of queries contact Customer Service (see section 2.3).

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see Fig. 8: AUTOCON® III 400 UH 400.

3.2.2 Symbols on the packaging

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. It can be downloaded from www.karlstorz.com.

3.3 Scope of supply*

- 1 AUTOCON® III 400 UH 400 or UH 401 (UH 401 incl. option BiVascularSafe)
- 1 Power cord 400 A
- 1 KARL STORZ USB-Stick 20 0402 82
- 1 Instruction manual with inspection protocol
- * also see page 106 (specifications).

3.4 Components required for operation

- Power cord
- Footswitch
- Neutral electrode for unipolar applications
- Connecting cable for neutral electrode or instrument
- Instrument (unipolar or bipolar)

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4 Preparation

4.1 Setting up the HF device

NOTE: HF equipment generates high-frequency electromagnetic fields that can interfere with the operation of sensitive electronic equipment and, for example, lead to image distortion. HF equipment, and particularly its cabling, should thus never be installed close to any equipment sensitive to electromagnetic interference. Parallel cable routing with cables from other devices should be avoided as a general rule, as HF cables emit HF emissions when activated, which may interfere with other equipment. Additional measures include:

- Connect the power connection of the AUTOCON® III 400 or the impaired device to another circuit.
- Increase the distance between the AUTOCON® III 400 and the impaired device.
- Lay instrument cables so that, as far as possible, they are not in the vicinity of other devices and their connecting cables.

WARNING: Risk of electric shock! Always connect the HF device to a mains power system with a protective conductor in order to prevent electric shock.

DANGER: Risk of burns to patients due to excessive leakage current! Locate the HF device outside the immediate vicinity of the patient, see section Ambient conditions, page 7.

HF devices may be used only in rooms used for medical purposes that meet the requirements of DIN VDE 0100-710 or IEC 60364-7-710.

If the HF device was previously stored or transported at temperatures below +10 °C or a relative humidity above 75 %, non-condensing, it will take approximately three hours to adjust to room temperature.

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Preparation

1. Observe the specified operating conditions, see section Operating conditions, page 21.
2. Place the HF device on one of the following platforms:
 - a table,
 - an equipment cart,
 - a console suspended from a ceiling support or wall-mounted brackets.
3. Place the HF device a sufficient distance away from other electronic equipment, see section EMC, page 221.
4. Position the HF device with the front of the device facing the patient and surgeon.
5. Do not place any other device on the HF device.
6. Do not place any other objects on or above the HF device.
7. If necessary, place the HF device on the Argon beamer.
8. Connect the potential equalization cable to the potential equalization connector. The device's ground line should be installed by a qualified electrician.
9. Connect the power cord.

4.2 Switching on the HF device

i Do not use the HF device if the display components are not working! For troubleshooting see section Recognizing and remedying errors, page 79.

1. Switch the HF device on using the power switch on the rear side of the device and then touch the permanently lit Standby button on the front panel.
- The HF device performs a self-test: All display and control elements light up.
2. Check all display and control elements for proper operation:
 - Power switch
 - Touch screen
 - Unipolar connection sockets
 - Bipolar connection sockets
 - Activation bar for unipolar and bipolar sockets

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Preparation

- ☛ The main screen appears, and the HF device is ready for use.
- ☛ The parameters of the most recently selected program appear on the display.
- ❗ If the device is switched off at the front panel, the Standby mode is activated.

4.3 Connecting instruments

▷ Before connecting instruments, ensure that the following conditions are met:

- Combinations of accessories not mentioned in the instruction manual may be used only if they are explicitly designed for the intended use. Always observe performance characteristics and safety requirements.
- The insulation of the accessories (e.g. HF cables and instruments) must be sufficient for the maximum peak output voltage (see IEC 60601-2-2 and IEC 60601-2-18).
- The HF device is designed for connection of HF cables measuring 3 m or 4 m in length.
- Do not use accessories with defective insulation.

4.3.1 Instruments for unipolar application

1. Plug the neutral electrode cable into the socket for the neutral electrode and choose the corresponding neutral electrode type, see chapter Selecting the neutral electrode, page 36.
 - ☛ The socket illumination goes dark.
 - ☛ The neutral electrode button changes from gray to the measurement color (green, yellow or red).
2. Connect the electrode handpiece to one of the two unipolar connection sockets.
 - or –
 - In the case of an accessory without a finger button: Connect the footswitch and unipolar connecting cable.
 - or –
 - Connect the unipolar cable to one of the two unipolar connection sockets for unipolar instruments.

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4.3.2 Instruments for bipolar application

1. Connect the bipolar cable to the instrument (e.g. forceps).
 2. Connect the bipolar cable to one of the two bipolar connection sockets.
 3. For bipolar use without AUTOSTART, connect the footswitch.
– or –
Select a mode with AUTOSTART function at the corresponding socket.
- ☛ Once the instrument is connected, the application starts after the configured reaction time.

i When introducing bipolar instruments into a metal trocar, the AUTOSTART function should be disabled beforehand!

4.3.3 Connecting the footswitch

In addition to the manual switch, a footswitch can be used to activate various operating modes.

- ▷ Connect the desired footswitch only during operation to one of the two connection sockets for footswitches.

☛ The HF device automatically detects the connected footswitch and indicates this on the front panel display, including the selected connection socket.

i One dual-pedal footswitch and one single-pedal footswitch can be connected. Footswitches without a changeover switch cannot be used (see table).
During the operation, in addition to the footswitch, only the connection to the Argon beamer with fiber-optic cables and, if relevant, the connection to the KARL STORZ OR1™ control NEO system may be plugged into the rear of the device.

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The following footswitches can be connected to the HF device:

Art. no.	Designation
UF 901	Single-pedal footswitch with changeover function switch
UF 902	Dual-pedal footswitch with changeover function switch

4.4 Test for proper functioning

4.4.1 Automatic test function

The HF device automatically performs cyclic testing during operation. If any errors occur, see section Recognizing and remedying errors, page 79.

4.4.2 Performing the test for proper functioning

Perform the following functional test before putting the device into service:

 The accessories must be designed for the specified maximum voltage.

1. Connect the neutral electrode and attach it securely to the patient's arm.
 The EASY neutral electrode indicator changes to green.
2. Remove the neutral electrode.
 The indicator changes to red, acoustic signals sound.
-  The neutral electrode used for this test may not later be used for an operation.
3. Connect a unipolar HF handpiece to a unipolar socket if there is a green EASY indicator and use the manual switch and footswitch to individually activate 'Cut' and 'Coag'.
4. Check the settings on the display.
5. Now change to the bipolar output and connect bipolar forceps.
6. Select a mode with AUTOSTART, grasp moist gauze with the forceps, and check the display.

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7. Now change to a mode without AUTOSTART and use the footswitch to activate the bipolar output. Check the settings and indicators in the bipolar section.

4.4.3 Actions in case of problems

Proceed as follows in case of functional problems:

1. Immediately disconnect the patient from the HF device.
2. Perform a technical inspection of the HF device.
3. Observe the relevant national regulations. For example, in Germany report incidents and near-accidents to the German Federal Institute for Drugs and Medical Devices in accordance with Section 3 of the German Ordinance on the Installation, Operation and Use of Medical Devices (MPBetreibV). Observe the provisions of the in-house reporting system in this regard.
4. Consult the Technical Support, see section Technical Support, page. 100.

 The HF device can be switched off at any time using the power switch  on the rear of the device.

4.5 Neutral electrode monitoring

 Always use the largest possible electrode when attaching a neutral electrode.

4.5.1 General information

 KARL STORZ recommends using split neutral electrodes, since only this type of electrode allows the HF device to detect detachment of the neutral electrode if this occurs.

Monitoring of the neutral electrode minimizes the risk of burns at the site where the neutral electrode is attached.

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Two types of neutral electrodes can be monitored:

- Split electrodes for infants (for use with reduced power)
- Split neutral electrodes

The type of neutral electrode and its contact quality are selected and/or shown in the Neutral electrodes menu.

Errors related to the neutral electrode and possible remedies are shown on the display, see section Recognizing and remedying errors, page 79 (see Fig. 9, neutral electrode modes).

4.5.2 EASY neutral electrode monitoring (EASY monitoring)

i The maximum output power of unipolar current types is reduced to 50 W when a children's electrode is selected.

The EASY monitoring function measures changes in the resistance between the patient and the high frequency surgical unit before and during HF activation. If necessary, it generates visual and audible alarms to request staff intervention. This requires using a split neutral electrode with appropriate contact areas and suitable transition resistance, attached to the patient according to the manufacturer's instructions. The EASY system does not monitor the currents through the individual contact surfaces of the split neutral electrode.

A KARL STORZ electrode with a surface area of at least 90 cm² must be used for the unipolar program 'Resection' and the 'Moderate Coagulation' mode.

If an error message is generated, the display changes from green via yellow to red, depending on the type of error.

! **NOTE:** Incorrect application of the neutral electrode!
Ensure compliance with the specifications for correct attachment of the neutral electrode with regard to size, adhesive properties and full-surface contact of the complete electrode.

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5 Operation

5.1 Connecting the power supply

- ✓ The line voltage must match the voltage specified on the type plate.
 - ✓ The potential equalization cable must be connected to the potential equalization connector (see section 4.1). The device's ground line should be installed by a qualified electrician.
1. Connect the power cord to the generator and the plug to the grounded wall socket.
 2. Use the power switch 31 on the rear to switch the device on (see Fig. 10).
 - The Standby button on the front panel is permanently illuminated.
 3. Switch the device on at the Standby button.
 - The device performs a self-test.
 - The bars of the activation displays light up as do the white socket surrounds 11, 12, 13, 17 and 18 (see page 17).
 - Full functionality of the loudspeaker is indicated by the start melody.
 - The individually configurable start screen appears if it has been set up.

5.2 Program overview

5.2.1 Display

The status bar is located at the top of the display (see Fig. 11). The settings of the four sockets are shown below the status bar.

Settings can be configured for each of the sockets.

1. The 'Effect' button is used to set the effect of electrosurgical cutting or coagulation.
2. Assign the footswitch using the 'Pedal' button.
3. The 'Mode' button allows the desired type of current to be selected.
4. Click on 'max. Watt' to adjust the maximum output power.

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5.2.2 Status bar

The program name and the buttons 'Favorite', 'Key lock', 'EASY', 'Information' and 'Menu' are arranged in the status bar (see Fig. 12).

In combination with an Argon beamer and selection of an Argon mode, the 'Argon' button is also additionally displayed (see Fig. 13).

5.3 Activating and deactivating sockets

▷ To activate a deactivated socket, plug a connecting cable into the connection socket.

– or –

Press the button 'activate socket' (see Fig. 14).

☛ An overview of the socket settings appears.

The overview is grayed out if no instrument is connected to the socket (see Fig. 15).

The socket illumination extinguishes and the socket selection window lights up when an instrument is plugged in.

▷ To hide sockets, click on the 'off' button next to the socket setting overview (see Fig. 16).

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5.4 Unlocking the screen

The device screen locks automatically. To unlock it, click on any control element and push the slider which then appears from left to right. An open lock then appears in the status bar (see Fig. 17, Unlocking the screen).

To deactivate the key lock or to change the time up to automatic screen lock, see section Menu 'System settings', page 67.

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5.5 Configuring output currents

- i** All selection windows are closed after 30 seconds without assumption of the modification.
- i** If selection windows are open, all buttons outside this field are inactive. In this case, activations are not possible.
- i** A change to the currently loaded program, e.g. by adjusting the output, is indicated by the information 'changed' below the program name.

5.5.1 Selecting the mode


1. To select the type of current for unipolar cutting, first select one of the two sockets on the left side (see Fig. 18).
2. Press the yellow icon under the 'Mode' button.
 - A selection screen appears for the available modes, and the rim of the associated socket starts flashing.
3. Select the desired mode by pressing the corresponding button.
 - or –
 - Deactivate the mode by pressing the button 'Cutting off (X)'.
4. Press the '?' button for more information on this selection.
5. For additional selection options press ◀/▶ inside the selection window.
6. Confirm your selection by pressing the 'OK' button.
 - The main screen will be displayed.
 - or –
 - Press the 'Back' button to return to the main screen without changing the selection.
- i** If a mode is changed within a socket, the set parameters, e.g. effect and max. Watt, remain the same for the respective mode. However, if, for example, the factory default setting of a mode is adapted and subsequently changed to a different mode and then back again, the user changes are not undone.

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

5.5.2 Specifying power limits

1. To select the maximum output for unipolar cutting, first select one of the two unipolar sockets and press the yellow icon under the 'max. Watt' button (see Fig. 19).
2. Use the '+' and '-' buttons to adjust the power level in single steps.
– or –
Use the slider to set the power level in steps of 10.
3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.
– or –
Press the 'Back' button to return to the main screen without changing the selection.

5.5.3 Selecting the effect

1. To select the effect for unipolar cutting, first select one of the two unipolar sockets and press the yellow icon under the 'Effect' button (see Fig. 20).
2. Use the '+' and '-' buttons to adjust the effect in individual steps.
– or –
Use the slider to set the effect.
3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.
– or –
Press the 'Back' button to return to the main screen without changing the selection.

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

22

5.5.4 Assigning the footswitch

i Handpieces and instruments with manual switches can be activated without a configuration setting.

A dual-pedal footswitch and/or single-pedal footswitch, each with a changeover switch, can be connected.

The changeover switch enables switching between pedal levels.

1. Call up the footswitch menu by pressing the button 'Pedal' (see Fig. 21).
2. Select the desired footswitch by pressing the corresponding button. For example, choose the active pedal level for cutting and coagulation for the upper left socket.
– or –
Deactivate the footswitch by pressing the button marked with an 'X'.

The edge of the selected button lights up green (see arrow in Fig. 21).

3. Confirm the selection by pressing the 'OK' button.
– or –
Press the 'Back' button to return to the main screen without changing the selection.

The socket is assigned to the active pedal level.

4. Pedal levels can be changed using the footswitches. Press the black changeover switch to change the socket.

The orange background of the changeover switch in the Pedal button indicates that the lower left-hand socket is activated (see Fig. 22).

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i If two footswitches are connected, either a single-pedal footswitch or a dual-pedal footswitch can be selected for coagulation (see Fig. 23).

The following pedal icons can be distinguished between:

Icon/ Button	Description	Icon/ Button	Description
	Dual-pedal footswitch CUT active		Dual-pedal footswitch COAG active
	Dual-pedal footswitch CUT inactive		Dual-pedal footswitch COAG inactive
	Dual-pedal footswitch not connected		Dual-pedal footswitch not connected
	Single-pedal footswitch COAG active		Single-pedal footswitch not connected
	Single-pedal footswitch COAG inactive		Deactivate footswitch CUT
	Deactivate footswitch coagulation		

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i The maximum power output of unipolar current types is reduced to 50 W when a children's electrode is selected.

2. Select the type of connected neutral electrode by pressing the corresponding symbol.
EASY: for monitoring split neutral electrodes (see Fig. 24)
BABY: for monitoring split neutral electrodes for infants (see Fig. 25)
MONO: to select a one-piece neutral electrode (see Fig. 26)

3. Press the '?' button for more information on this selection.

4. Confirm your selection by pressing the 'OK' button.
– or –
Press the 'Back' button to return to the main screen without changing the selection.

The selected type of neutral electrode in connection with a color indicator for the contact quality is shown in the status bar.

i When using the 'EASY' and 'BABY' mode, no unsplit electrodes are accepted. Using the 'MONO' mode, no split electrodes are accepted.

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According to the contact quality, the icons are shown in color for neutral electrodes:

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Icon/ Button	Description
	Split neutral electrode contact quality OK
	Split neutral electrode contact quality not optimum
	Split neutral electrode contact quality insufficient
	Split neutral electrode not connected
	Split baby neutral electrode contact quality OK
	Split baby neutral electrode contact quality not optimum
	Split baby neutral electrode contact quality insufficient
	Split baby neutral electrode not connected
	Non split neutral electrode contact quality OK
	Non split neutral electrode not detected or contact quality insufficient
	Non split neutral electrode not connected
	Display contact quality

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5.5.6 Coding system

The coding system enables automatic instrument identification. The coding system recognizes the connected coded KARL STORZ instrument and automatically selects the preferred parameters.

1. Insert the coded instrument in an AUTOCON® III 400 socket.

☛ The instrument data is read (see Fig. 27).

☛ A description of the instrument appears:

- Instrument name
- Recognized socket
- Article number
- Batch number
- Remaining usages

☛ The parameters are transferred automatically after 5 seconds and shown in the main screen.

• If the coded instrument is connected to a socket without preset parameters, with the coded system the ideal settings are loaded for the coded KARL STORZ instrument.

• If the coded instrument is connected to a socket with preset parameters, a plausibility check is performed. If the preset values for the coded instrument are within the permitted range, these are not overwritten. The coded instrument can now be used with the preset parameters.

– or –

The selection is saved by pressing the 'OK' button.

To return to the main screen without changing the selection, press the 'back' button.

☛ The coded instrument can now be used.

☛ The permitted parameters for the coded KARL STORZ instrument can still be selected, all other modes are gray.

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5.6 Mode overview

An overview of the current types that can be executed with the HF device is shown below.

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Cutting mode symbol	Designation
	Gastro Knife 2
	Gastro Knife 3
	Argon*

Coagulation mode symbol	Designation
	Laparoscopy
	Moderate
	Forced coag
	Resection
	Spray
	Forced mixed
	Forced cutting
	Gastro Coag
	Argon flexible*

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Coagulation mode symbol	Designation
	Argon flex. pulse*
	Argon open*
	Cardiac Thorax
	Cardiac Mammary
	SimCoag

i * These modes can only be activated in combination with an Argon beamer. Our service team will be happy to provide you with information on the available combination options.

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5.6.2 Bipolar modes

Cutting mode symbol	Designation
	Bip. cutting
	Bipolar scissors
	Bip. resection
	Bip. vaporization

Coagulation mode symbol	Designation
	RoBi®
	Laparoscopy
	Standard
	Bip. resection
	Bip. vaporization
	BiVascularSafe*
	BiVascularSafe Open Surgery*

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Coagulation mode symbol	Designation
	Standard AUTO
	Bipolar scissors
	Micro
	Forced

i * These modes are available for the device versions UH 401 and UH 401U.

The information and data regarding settings, application points, application duration and instrument use are based on clinical practice. However, these are only basic guidelines which must be tested for suitability by the operator. Depending on individual conditions, it may be necessary to deviate from the provided data. Medical practice is continuously evolving as a result of R&D and clinical experience. This may also make deviations from the provided data necessary.

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5.7 Unipolar cutting modes

5.7.1 Laparoscopy

This mode is used in laparoscopy and arthroscopy for unipolar cutting.

Application areas

Laparoscopy, arthroscopy

Suitable instruments

- Arthroscopy electrodes
- Laparoscopy electrodes

5.7.2 Standard

In this mode a high-performance HF current with a low crest factor is used for cutting biological tissue.

ARC CONTROL adjusts the power output to the minimum required level in response to variations in tissue type and changes in the cutting area or speed.

Application areas

Cutting tissue with low electrical resistance, such as muscle tissue or vascular tissue.

Cutting or preparing fine structures

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- Electrodes for arthroscopy
- Electrodes for laparoscopy

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- Resektoskop (unipolar)
- Resektionsschlinge
- Vaporisations-Elektrode

5.7.3 Micro

This mode is used for electrosurgical cutting using micro-electrodes. It enables extremely fine control of the power level and precise work.

Application areas

Pediatric surgery, neurosurgery,
plastic surgery

Suitable instruments

- Micro needle electrodes

5.7.4 Resection

This mode is used in gynecology and urology. Arc control generates the cutting effect with simultaneously minimized output power. ARC Control facilitates direct cutting and prevents electrode adhesion.

- i** Use non-conductive irrigation fluids.

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- Gynäkologische Schlingen

5.7.5 Resection C-Cut®

This mode is used for electrosurgical cutting using unipolar cutting instruments and electrodes.

The cut features a particularly pronounced level of hemostasis.

Different hemostasis effects can be set.

ⓘ Use non-conductive irrigation fluids.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), vaporization of prostate tissue (TUR-VAP).

Suitable instruments

- Resectoscope (unipolar)
- Resection loop

5.7.6 SupraLOOP

This mode is used in gynecology for laparoscopic removal of the uterus. Removal can be achieved by emitting a unipolar cutting current and, at the same time, tensioning the loop.

Application areas

Gynecological, laparoscopic removal of the uterus.

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- Polypektomieschlingen

5.7.7 Dry

This mode is used for unipolar dry cutting. A large, controlled arc is generated, which allows significantly deeper coagulation to be obtained.

Application areas

Cardiac surgery and blood coagulation in retracting blood vessels in the sternum region.

Suitable instruments

- Knife electrodes

5.7.8 Gastro Loop 1

This mode is used in gastroenterology. Polypectomy snares are used for cutting and coagulation. Arc control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a relatively slow pulse rate of 1 cutting pulse per second, this mode is suitable for especially cautious work.

Application areas

Polyp removal using polypectomy snares and flexible endoscopy.

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- Polypektomieschlingen

5.7.9 Gastro Loop 2

This mode is used in gastroenterology.

Polypectomy snares are used for cutting and coagulation. Arc control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic pulse rate of 1.5 cutting pulses per second this mode is suited to experienced users.

Application areas

Polyp removal using polypectomy snares and flexible endoscopy, with dynamic pulse rate for experienced users.

Suitable instruments

- Polypectomy snares

5.7.10 Gastro Loop 3

This mode is used in gastroenterology.

Polypectomy snares are used for cutting and coagulation. Arc control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic and fast pulse rate of 2.2 cutting pulses per second, this mode is suitable for advanced users.


Application areas

Polyp removal using polypectomy snares and flexible endoscopy, with dynamic fast pulse rate for advanced users.

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5.7.11 Gastro Knife 1

This mode is used in gastroenterology. Instruments for papillotomy and endoscopic resections are used for cutting and coagulation. Arc control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a relatively slow pulse rate of 1.3 cutting pulses per second, this mode is suitable for especially cautious work.

Application areas

Papilla incision using a papillotome and flexible endoscopy, resection with needle knives; slow pulse rate for cautious work.

Suitable instruments

- Papillotomes
- Needle knives

5.7.12 Gastro Knife 2

This mode is used in gastroenterology. Instruments for papillotomy and endoscopic resections are used for cutting and coagulation. Arc control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic pulse rate of 1.8 cutting pulses per second, this mode is suitable for experienced users.

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- Papillotome
- Nadelmesser

Application areas

Papilla incision using a papillotome and flexible endoscopy, resection with needle knives; dynamic pulse rate for experienced users.

Suitable instruments

- Papillotomes
- Needle knives

5.7.13 Gastro Knife 3

This mode is used in gastroenterology. Instruments for papillotomy and endoscopic resections are used for cutting and coagulation. Arc control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic and fast pulse rate of 2.2 cutting pulses per second, this mode is suitable for advanced users.

Application areas

Papilla incision using a papillotome and flexible endoscopy, resection with needle knives; dynamic fast pulse rate for advanced users.

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5.7.14 Argon

This mode is used to perform open surgical interventions in combination* with an Argon beamer for argon-assisted cutting. With suitable instruments connected, argon-assisted cutting can be performed using rigid electrodes.

Application areas

Visceral surgery

Suitable instruments

- Rigid argon electrodes
- Argon handpiece

 * Our service team will be happy to provide you with information on the available combination options.

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5.8 Unipolar coagulation modes

5.8.1 Laparoscopy

This mode is used in laparoscopy and arthroscopy for unipolar coagulation.

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- Kugelelektroden
- Messerelektroden
- Spatelektroden

5.8.2 Moderate

This mode is used with contact coagulation to stop hemorrhagic oozing, for hemostasis of relatively large tissue areas, and for small-area coagulation. Tissue carbonization is prevented and electrode adhesion to the tissue is strongly reduced. Greater coagulation depth is achieved in this mode than in other coagulation modes. The degree of surface scabbing can be controlled by adjusting the 'Effect' setting in the range of 1 to 3.

Application areas

Coagulation with high penetration depth; low electrode adhesion to tissue.

Suitable instruments

- Electrodes with large contact areas, such as ball electrodes

5.8.3 Forced coag

This mode is used for contact coagulation with low tissue penetration, preferably using fine electrodes and electrodes with small contact areas. It achieves a high degree of coagulation with low cutting tendency.

Application areas

Fast coagulation with small penetration depth

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- Messerelektroden
- Spatelektroden
- Nadelelektroden

5.8.6 Forced mixed

This mode is used for contact coagulation with low tissue penetration, preferably using fine electrodes and electrodes with small contact areas. It achieves a high degree of coagulation with moderate cutting tendency.

Application areas

Fast coagulation with small penetration depth and moderate cutting tendency.

Suitable instruments

- Knife electrodes
- Spatula electrodes
- Insulated unipolar forceps

5.8.7 Forced cutting

This mode is used for contact coagulation with low tissue penetration, preferably using fine electrodes and electrodes with small contact areas. It achieves good hemostasis with very good cutting tendency.

Application areas

Fast coagulation with small penetration depth and very good cutting tendency.

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5.8.10 Argon flex. pulse

This mode is used for argon-assisted electrosurgery in combination* with an Argon beamer. This is the current type Spray. For argon-assisted coagulation, flexible probes are used in combination with endoscopes. The pulse frequency changes with the effect setting. The higher the effect level, the faster the pulse sequence.

Application areas

Gastroenterology, homogeneous surface coagulation

Suitable instruments

- Flexible argon probes

5.8.11 Argon open

This mode is used to perform open surgical interventions in combination* with an Argon beamer for argon-assisted electrocoagulation.

This is the current type Spray. With suitable instruments connected, argon-assisted coagulation can be performed using rigid electrodes.

Application areas

Visceral surgery

Suitable instruments

- Rigid argon electrodes
- Argon handpiece

* Our service team will be happy to provide you with information on the available combination options.

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- Kugelelektroden
- Messerelektroden
- Spatelektroden

5.8.14 SimCoag

This mode is used for simultaneous preparation. Two unipolar sockets can be activated at the same time to enable the simultaneous use of two manually switched instruments. Both handpieces can be switched on and off independently.

This mode is used with non-contact surface coagulation using an arc.

i The output power setting applies to both output sockets, and the power is distributed according to the tissue structure.

! **WARNING:** The output power of all active electrodes can change during use.

Application areas

Simultaneous coagulation and preparation, e.g. for cardiac or mammary surgery

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
5.9.3 Bip. resection

This bipolar mode is used in gynecology and in urology for resection with loop electrodes under conductive irrigant solution (saline). Arc control generates the cutting effect with simultaneously minimized output power. ARC Control facilitates direct cutting and prevents electrode adhesion.

i Make sure that NaCl is used as an irrigation medium.

Secure a continuous irrigation during the application.

Always use conductive lubricants to avoid damages of the urethra.

Avoid continuous activations.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT).

Suitable instruments

- Resectoscope (bipolar)
- Resection loop

i Optimum results are possible only when using the coded KARL STORZ resection cable.

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
5.9.4 Bip. vaporization

This bipolar mode is used in gynecology and urology with vaporization. Upon tissue contact the arc is ignited directly resulting in rapid tissue vaporization with minimal heat emission to the surrounding area.

- i** Make sure that NaCl is used as an irrigation medium.
Secure a continuous irrigation during the application.
Always use conductive lubricants to avoid damages of the urethra.
Avoid continuous activations.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT), vaporization of prostate tissue (TUR-VAP).

Suitable instruments

- Resectoscope (bipolar)
- Vaporization electrode

Bipolar scissors should only be operated with the current type bipolar scissors for cutting or coagulation.

5.10 Bipolar coagulation modes

5.10.1 RoBi®

This mode designed specially for RoBi® is used for coagulation in combination with bipolar laparoscopic instruments.

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- i** Make sure that NaCl is used as an irrigation medium.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT)

Suitable instruments

- Resectoscope
- Resection loop
- Vaporization electrode

- i** Make sure that the instrument has contact with the tissue while activating bipolar coagulation to avoid an unintended heating of the irrigation fluid.

5.10.5 Bip. vaporization

This mode is used for bipolar hemostasis in gynecology and in urology with vaporization.

- i** Make sure that NaCl is used as an irrigation medium.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT), vaporization of prostate tissue (TUR-VAP)

Suitable instruments

- Resectoscope
- Vaporization electrode

- i** Make sure that the instrument has contact with the tissue while activating bipolar coagulation to avoid an unintended heating of the irrigation fluid.

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**5.10.6 BiVascularSafe
(only UH 401/UH 401U)**

This mode is used for the permanent sealing of veins, arteries and tissue bundles. Conventional instruments cannot be used in this mode. All parameters are preconfigured and automatically regulated.

Application areas

Vessel sealing open and laparoscopic


Suitable instruments

- Sealing instruments*

Suitable electrodes*

 *Please contact your KARL STORZ authorized dealer for further advice and information.

**5.10.7 BiVascularSafe Open Surgery
(only UH 401/UH 401U)**

This mode is used for the permanent sealing of veins, arteries and tissue bundles with open surgical sealing instruments. Conventional instruments cannot be used in this mode. All parameters are preconfigured and automatically regulated.

Application areas

Vessel sealing in open surgery

Suitable instrument

- Open surgery sealing instruments*

 *Please contact your KARL STORZ authorized dealer for further advice and information.

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5.10.8 Standard AUTO

This mode is used for arcless contact coagulation with forceps. Activation starts automatically on contact with tissue. The adjustable delay time can be set under MENU – SYSTEM SETTINGS – AUTOSTART DELAY (page 67).

- i** Setting the AUTOSTART mode can result in unintentional coagulations, e.g. when bipolar forceps are used for gripping while the AUTOSTART mode is on.
- i** When introducing bipolar instruments into a metal trocar, the AUTOSTART function should be disabled beforehand!

Application areas

Bipolar coagulation with AUTOSTART

Suitable instruments

- Bipolar forceps

5.10.9 Bipolar scissors

This mode is used with bipolar scissors. It can be used for coagulation before or during mechanical cutting, point coagulation and surface coagulation.

Application areas

Preparation, coagulation and cutting of tissue

Suitable instruments

- Bipolar scissors

- i** Bipolar scissors should only be operated with the current type bipolar scissors for cutting or coagulation.

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- bipolare Pinzetten

5.10.10 Micro

This mode is used for arcless contact coagulation with micro forceps. It enables extremely fine control of power output down to 0.1 W and precise work for tightly restricted bipolar contact coagulation.

Application areas

Bipolar coagulation in pediatric surgery, neurosurgery, plastic surgery, etc.

Suitable instruments

- Bipolar forceps
- Micro forceps

5.10.11 Forced

This mode is used for forced coagulation with forceps.

i Given the higher voltage, we do not recommend the use of this current type for minimally invasive surgery. Note the permissible instrument voltage.

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5.11 Menu

i Basic settings can be changed, e.g., operational language, audio, display and saving options.

5.11.1 Overview

The functions shown in the figure on the left (Fig. 28) are available.

Selecting a menu

▷ Press the corresponding selection button to open the menu.

Exiting a menu

▷ Press the 'back' button to return to the main screen.

5.11.2 Menu 'System settings'

In the 'System settings' menu the settings shown in the figure on the left (Fig. 29) can be changed.

Autostart delay

The time can be set under Autostart delay after which the HF generator is automatically activated impedance-dependent following tissue contact (see Fig. 30).

Incremental setting

▷ Press the buttons '+' and '-'.

Fast setting

▷ Move the slider to and fro.

Language selection

The following languages can be selected under 'Language selection':

German, English, Spanish, French, Portuguese, Italian, Russian, Polish*, Turkish*, Japanese*, Korean*, Thai*, Indonesian*, Chinese*, Czech*, Arabic*, Hungarian*, Danish*, Finnish*, Vietnamese*, Swedish*, Norwegian*, Dutch*, Bulgarian*, Serbian*, Lithuanian*, Romanian*, Slovakian*
* at a later date

Key lock

In order to switch off the automatic screen lock or to set its duration press 'Key lock'. The duration can be set between 30 seconds and 5 minutes.

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OR system selection

The connection to the KARL STORZ OR system OR1™ can be created under the menu item 'OR system selection' (see Fig. 31). Upon delivery the KARL STORZ OR system OR1™ is in Standby mode. This is evident by the symbol .

To create a connection to your OR1™ a LAN cable specified by KARL STORZ must be used. You can get more detailed information on this from our service team. If the OR system OR1™ is in Standby mode, a connection is established automatically when the connecting cable is inserted.

 A detailed description can be found in the OR1™ manual SCB control NEO art. no. 96206572 DER.

Deactivating the OR system (operating room without OR1™ system):

1. Click on 'off'.
2. Press OK to confirm.

 The OR1™ symbol is hidden. An automatic connection upon connecting the OR1™ system is no longer possible.

Returning the OR system (after deactivation) to Standby mode

1. Select 'OR1, Karl Storz'.
2. Press OK to confirm.

 The OR1™ symbol is displayed .

OR1™ configuration (only for trained service personnel)

Special network configurations are possible by clicking on the symbol  (see Fig. 31).

 A detailed description can be found in the OR1™ manual SCB control NEO art. no. 96206572 DER.

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5.11.3 Menu 'Volume'

Use the 'Volume' menu to set the volume of the individual acoustic signals (see Fig. 32).

Incremental setting

▷ Press the buttons '+' and '-'.

Fast setting

▷ Move the slider to and fro.

i The volume of the activation signal should be increased as necessary for use in relatively noisy surroundings. The alarm tones have a minimum volume and limited changeability.

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Mode	Category	Frequency (Hz)	Signal type
Unipolar Cut	Activation tones	635	Continuous sound
Unipolar Coag	Activation tones	475	Continuous sound
Bipolar Cut	Activation tones	565	Continuous sound
Bipolar Coag	Activation tones	505	Continuous sound
Sim Coag	Activation tones	755	Continuous sound
BiVascularSafe end	Activation tones	-	Alternating sound
Footswitch changeover	Activation tones	-	Signal tone
ZAP mode	Activation tones	-	Signal tone
Fault	Alarm tones	-	Signal tone
Warning	Alarm tones	-	Signal tone
Note	Alarm tones	-	Signal tone

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5.11.4 Menu 'Service'

Service contact data can be called up in the 'Service' menu (see Fig. 33).

After entering a password, you can make further settings, e.g. resetting to the delivery state. To access the service level ('Service Tools') (see Fig. 34), enter the password 001224.

Saving device settings

With the function '1 Backup device' device settings can be saved on the KARL STORZ USB-Stick 20 0402 82. This includes all saved programs and system settings.

Installing device settings

Use '2 Restore device' to transfer saved device settings on the KARL STORZ USB-Stick to an AUTOCON® III 400.

Adding a start screen

Select the function '3 Add logo' to display a user-defined logo when starting up.

Deleting a start screen

Select the function '4 Remove logo' to remove a user-defined logo.

Program settings

Select the function '5 Program settings' to go to the setting level of the program/procedure configuration (see Fig. 35). In this menu restrictions ('Program restrictions') for the program handling and visibility of the procedures defined by the manufacturer can be set. Upon delivery, programs can be saved by the user (cannot be deleted or overwritten) and all KARL STORZ procedures are selected.

Saving logfiles

Select the function '6 Save logfiles' to save all logfiles onto the KARL STORZ USB-Stick.

Resetting to default setting

Select the function '7 Reset to default' to reset all settings and programs to the default setting.

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5.11.5 Menu 'System information'

In the 'System information' menu various system parameters such as version, serial number, TSI date for AUTOCON® III 400 and, if applicable, Argon-Beamer are displayed, as well as options (see Fig. 36).

For more details see the section Safety inspection, page 97.

5.11.6 Menu 'Programs'

1. In the 'Programs' menu you can select programs from a list and assign to or delete from favorites.
– or –
Click on the current program name in the main screen to access this menu quickly.
 2. Press the relevant program name to select the program.
 3. To navigate vertically in the program list, press the corresponding button ▲/▼. The programs are always sorted in alphabetical order.
 4. Assignment to favorites is possible using the 'star buttons' at the bottom of the screen. Add programs to favorites using the 'green arrow' star button ➤ (see Fig. 37 or 38 – left column) and remove these using the 'red arrow' star button ➤ (see Fig. 39). The button ➤ appears in place of the ➤ button if a program is marked in the left-hand column.
 5. Sort your favorites using the ➤/➤ star buttons.
 6. Load the selected program by pressing 'OK'.
– or –
Return to the main screen by pressing 'back'.
- ❗ A maximum of 200 favorites can be selected.

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Saved programs can be deleted in the 'Programs' menu (if the function is enabled under 'Program settings' in section 5.11.4).

1. To this end, select the program to be deleted in the program list by clicking on the program name (see Fig. 36 or 37). To navigate vertically in the program list, press the corresponding button ▲/▼.

2. To delete the selected program permanently, select the 'recycle bin' button.

☛ The selected program is deleted upon confirmation of a prompt.

5.11.7 Menu 'Favorites'

The defined favorites can be selected in the 'Favorites' menu. A fast selection of the favorites is possible using the star button in the main screen.

Access the next page of the favorites list using the buttons ◀/▶ in the lower part of the screen (see Fig. 40).

1. Confirm with 'OK' to accept the selection.
2. Press the 'back' button to return to the main screen.

5.11.8 Menu 'Save program'

In the 'Save program' menu (see overview under 5.11.1) the current setting can be saved under the same or another program name (see Fig. 41). This function is enabled as standard (upon delivery); see 'Program settings' in section 5.11.4.

Program names can be created using the keyboard. Symbols, capital or small letters or numbers are selection options.

Program names spread over two lines can be assigned using the 'Enter' button.

1. Confirm with 'OK' to accept the selection.
2. Press the 'back' button to return to the main screen.

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
5.11.9 Socket settings

In the 'Socket configuration' menu the settings shown in the left-hand figure (Fig. 42) can be changed.

Expanding sockets

1. Call up 'Expand socket' in the 'Socket configuration' menu (see Fig. 43).
2. Click on  next to the two connectors to change the current setting (to divide the second bipolar socket on the bottom of the device into 2; see page 17, number ⑩).

With active socket expansion the  on the left of the two connectors changes to .

3. Confirm with 'OK' to accept the selection.
 4. Press the 'back' button to return to the main screen.
- You can now use the bipolar socket on the bottom of the device as two sockets.

ZAP mode

With the ZAP mode you can switch between two preset current forms for the same instrument.

1. Call up the ZAP mode in the 'Socket configuration' menu (see Fig. 44).
2. Activate or deactivate the changeover for the individual sockets by pressing the 'tick' next to the relevant socket (see Fig. 45).

If the ZAP mode is active, the 'tick' is displayed filled in.

Confirm with 'OK' to accept the selection.

Press the 'back' button to return to the main screen.

In the main screen a second level of the socket settings is shown for the active ZAP mode (see Fig. 46).

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You can switch in the main screen using either the handpiece or foot pedal.

For manual switching in the main screen tap on the ZAP mode symbol next to the socket settings (see Fig. 45).

On changing the ZAP mode level, the color of the symbol changes from *black* to *white* (see Fig. 46).

Another way of changing the ZAP mode level is by actuating both buttons on the handpiece.

You can also change levels via the black changeover switch on the footswitch.

▷ To do this, select the ZAP mode under the 'Pedal' icon (see Fig. 47).

✦ The ZAP mode symbol now appears in the main screen below the 'Pedal' icon (see Fig. 48).

▷ To change the ZAP mode levels, actuate the black switch on the footswitch.

5.11.10 Menu 'System messages'

The 'System messages' menu can be accessed in the main screen via the 'Information' symbol.

Here the errors saved in the HF device which have occurred since switching on the HF device are called up (see Fig. 49).

These messages are not saved when switching off the HF device.

Calling up saved system messages:

1. Select a system message.
2. The selected system message is displayed again with '?'.
3. Press 'OK' to return to the overview.

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Argon flexibel:
COAG Flussrate: 0,4 l/min

5.11.11 Menu 'Argon'

In case of the successful recognition of an Argon beamer device, this menu is selectable in the status bar.

The 'Argon' menu enables the selection of argon flow rates for cutting and coagulation, as well as the selection of argon bottles and the display of filling levels.

The Argon menu can be selected quickly via the 'Argon' button in the main screen (see Fig. 50).

1. Use the '+' and '-' buttons to adjust the argon flow rate for cutting (Cut) or coagulation (Coag).
– or –
Use the slider to set the flow rate.
2. In case of two connected argon bottles use the buttons for output '1' or '2' to select the desired gas source.

Pressure reducers with electronic pressure sensor enable the display of exact filling levels of argon bottles.

3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.
– or –
Press the 'Back' button to return to the main screen without changing the selection.

i Default settings for argon flow rates according to the different modes are selected automatically:

Argon open:
CUT flow rate: 4.0 l/min
COAG flow rate: 3.0 l/min

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


5.11.12 Basic programs (procedures)

Procedures are programs in which parameters such as voltage and output are determined for specialist medical interventions. You can either create your own procedures/programs (see section 5.11.8) or use procedures defined by KARL STORZ. These procedures are assigned to different disciplines, e.g. urology, laparoscopy and start with the number 0100. All specialist KARL STORZ procedures/programs and a standard program are enabled (visible) as standard (upon delivery). These programs and newly created ones cannot be changed or deleted (unless the corresponding setting in the menu Service – Program Restrictions, section 5.11.4, has been changed). In this menu you can display/ hide the specialist procedures to suit your needs. For selection see section 5.11.6 Menu 'Programs'.

i The procedure name is displayed in the status bar (see section 5.2.2).

The predefined procedures/programs set in the factory for the AUTOCON® III 400 are divided into the following user groups:

- Specialization urology (0100-0105; see Fig. 51 to 56)
- Specialization gynecology (0200-0208; see Fig. 57 to 64)
- Specialization laparoscopy (0300-0303; see Fig. 65 to 68); Procedure 0303: BiVascularSafe only possible with models UH 401 and UH 401U
- Specialization gastroenterology (0400-0401; see Fig. 69 to 70)
- Application open surgery (0500-0501; see Fig. 71 to 72) Procedure 0501: BiVascularSafe only possible with models UH 401 and UH 401U
- Specialization pediatrics (0600; see Fig. 73)

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Operation – Overview of procedures












- Specialization orthopedics and sports medicine (0700; see Fig. 74)
- Specialization ENT (0800; see Fig. 75)
- Specialization visceral surgery (0900; see Fig. 76); BiVascularSafe only possible with device models UH401 and UH401U

i An overview (screenshots) of the predefined procedures can be found on this page and the following page.

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








5. 11. 13 Switching off the HF device (decommissioning)

1. Switch off the device at the Standby button on the front panel.
2. Switch off the device at the power switch ❶ on the rear (see Fig. 10).
3. Remove connected accessories from the device and, if necessary, reprocess.

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❶ When the power switch on the rear is turned off, the Standby button (surrounding light) remains white for a few seconds. Once the surrounding light has gone off, the device is voltage-free.

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- System error
- EASY monitoring error

6.1 System information

System information is displayed as a message on the screen (see Fig. 77).

For further information on the cause and how to remedy this message, press the 'i' button in the status bar.

System information can be divided into three categories:

- Notices (gray display)
- Warnings (orange display)
- Faults (red display)

Notices are shown for five seconds on the screen. Warnings and faults are shown for ten seconds.

While a fault is present, activations are prohibited.

This information can be called up again under the now orange-framed 'i' in the status bar.

These messages are also available in the menu system information and will be deleted when switching off the device.

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The following table describes the cause of the error and the appropriate corrective action.

Header	Message text
Confirmation of AUTOSTART mode	You have selected an AUTOSTART mode. Setting the AUTOSTART mode can result in unintentional coagulations, e.g. when bipolar forceps are used for gripping while the AUTOSTART mode is on.
AUTOSTART Fault	The instrument is in contact with tissue. AUTOSTART cannot be selected when the instrument is in contact with tissue. Open up the instrument.
TSI Warning	The annual Technical Safety Inspection (TSI) is due.
Neutral Electrode Fault	No neutral electrode connected. No neutral electrode is connected. Connect a neutral electrode.
Neutral Electrode Fault	Wrong neutral electrode connected. The selection does not match the neutral electrode connected. Connect the neutral electrode which matches the selected mode, or change the mode to match the neutral electrode.
Neutral Electrode Fault	Poor contact with the patient. The resistance between the neutral electrode and the patient is too high. Establish better contact of the neutral electrode.
Neutral Electrode Warning	Poor contact with the patient. The contact resistance between the neutral electrode and the patient is getting worse. Establish better contact of the neutral electrode.
Neutral Electrode Fault	No cable for neutral electrode connected. No neutral electrode cable is connected. Connect a neutral electrode.
Neutral Electrode Notice	No cable for neutral electrode connected. The cable for neutral electrode has been removed. Unipolar activation is not possible.
Mode Fault	No mode selected. No mode was selected for this type of activation. Select the desired mode or change the footswitch assignment.
Mode Fault	This mode is not allowed for baby neutral electrodes. Use split neutral electrodes with a large surface and the suitable mode for this.
Mode Fault	This mode is not allowed for this socket. The current mode remains active. Choose another socket for this mode.
Footswitch Fault	No compatible footswitch. The connected footswitch is not compatible with this device. Connect a compatible footswitch with a changeover switch.
Footswitch Fault	Footswitch not assigned to a socket. The footswitch has not been assigned to a socket yet. Assign a socket to the footswitch using the 'Pedal' button.
Footswitch Fault	Error on footswitch connection. Check the footswitch. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.

Header	Message text
Fingerswitch Fault	Fault on fingerswitch connection. Check the handpiece and the connection cable. Please replace them if damaged. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.
Temperature Warning	The temperature of the device is higher than normal. The temperature of the device is elevated. This leads to a reduction of the maximum power.
Limitation of Continuous Activation	The maximum activation time has been exceeded. Please only activate the generator in short intervals, in order to avoid harming the patient and damaging the connected instruments or the generator.
Activation Fault	While switching on the device, there is an activation by footswitch, fingerswitch or AUTOSTART. Check the handpieces or footswitches for malfunctions. Disconnect the handpieces / footswitches from the device. If the error persists, please contact the Technical Support. Contact: MENU – SERVICE.
Activation Fault	There is an activation while connecting the footswitch or fingerswitch. Check the handpieces or footswitches for malfunctions. Disconnect the handpieces / footswitches from the device. If the error persists, please contact the Technical Support. Contact: MENU – SERVICE.
Activation Fault	There is no instrument connected on the activated socket. Connect an instrument on the designated socket.
Activation Warning	The mode for safety inspections is active. Activation is not possible. Quit this mode before activating again.
Bipolar Resection Warning	Use a coded KARL STORZ resection cable. Make sure that NaCl is used as an irrigation medium. Secure a continuous irrigation during the application. Always use conductive lubricants to avoid damages of the urethra. Avoid continuous activations.
GastroCut Warning	Polypectomy snare not in contact with tissue, or check connection cable at snare or generator. Please apply the snare and reactivate. First of all establish contact between tissue and polypectomy snare, or check the connection cable at the snare or the generator. Then activate with the yellow foot pedal.
BiVascularSafe Notice	Sealing has been completed.
BiVascularSafe Warning	There is a short-circuit in the area of the sealing instrument. Check the instrument for foreign bodies and whether there is contact with other objects. The sealing point must be free from foreign bodies.

Header	Message text
BIVascularSafe Warning	The sealing instrument is not in contact with tissue. Grasp the tissue and activate BIVascularSafe again. Check the connection to the generator. To seal the vascular and tissue bundles, it is necessary to grasp tissue and close the sealing instrument prior to activation.
Argon-Beamer Fault	Please connect the argon unit to the generator and switch it on. The argon device is connected to the generator with fiber optic cables. The active argon unit is automatically activated by the generator when switching to an argon mode.
Internal Error Argon- Beamer 5100	Argon-Beamer is inoperable. Connect an operational Argon unit to the generator. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.
Argon-Beamer Fault	Please check if the argon bottles are connected and open. Empty bottles should be replaced. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Fault	The argon inlet pressure is too high. Max. inlet pressure: <4.5 bar Connect a source of argon gas in the appropriate pressure range. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button.
Argon-Beamer Fault	The argon inlet pressure has exceeded the permissible limits. Inlet pressure range: 2 – 4.5 bar Connect a source of argon gas in the appropriate pressure range. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button.
Argon-Beamer Warning	Mixed operation of argon bottles with and without an electric bottle pressure gauge is not recommended. Connect two identical pressure reducers.
Argon-Beamer Warning	Please check if the instrument is free of adhesions, and purge it with argon. If repeated purging does not solve the problem, the instrument and cable must be replaced.
Argon-Beamer Fault	Please check if the argon bottles are connected and open. Empty bottles should be replaced. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Warning	The filling level of the argon bottle is low. Please make sure that a replacement is available. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Fault	The argon bottle is empty. Connect a replacement bottle to enable activation. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Notice	One argon bottle is empty. The unit has shifted automatically to the replacement bottle. Please make sure that a replacement is available.

Header	Message text
Argon-Beamer Fault	Please check if the argon bottles are connected and open. Empty bottles should be replaced. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button.
TSI Argon-Beamer Warning	The annual Technical Safety Inspection (TSI) for the Argon-Beamer is due.
Coded System Notice	The lifetime of the instrument is ending soon. Please order a replacement in good time. Any use of the instrument beyond its lifetime is not covered by warranty. Please contact your dealer in good time to purchase a new instrument.
Coded System Warning	The maximum lifetime of the instrument has been reached. Any further use is not covered by warranty. The maximum service lifetime of the instruments must not be exceeded, in order to guarantee safe usage. Any further use is at the user's risk.
Coded System Warning	A software update is necessary to use the Coded System with this instrument. Only carry out manual settings at this instrument. Please contact the Technical Support. Contact: MENU – SERVICE.
Coded System Warning	The preferred parameters for this coded instrument cannot be loaded. Only carry out manual settings at this instrument. Please contact the Technical Support. Contact: MENU – SERVICE.
Internal Error xxxx (e.g. with xxxx=4183)	If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.

With internal errors a number is shown in the header.

Please quote this number to the Technical Support.

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6.2 Error indications for EASY monitoring

Error indications are displayed in three stages (green, yellow and red) when problems occur.

When working with a split neutral electrode, the following errors may occur:

EASY monitoring	Cause	Indication	Remedy
Lights up yellow	Significant increase in resistance Depending on the indication, there may be heating under the neutral electrode	-	Stopping the application is not necessary. Check the positioning of the neutral electrode.
Switches from green to continuous red	A significant problem occurs when the unipolar current is activated.	An acoustic signal sounds. An error message appears on the display.	Check the neutral electrode and neutral electrode cable, see section EASY monitoring, page 28. ▷ Check the neutral electrode cable for proper connection and external damage.
	Loosened electrode	An acoustic signal sounds. A warning message appears on the display.	▷ Correct the positioning of the neutral electrode. If the error persists, replace the neutral electrode.

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! WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

! CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

! CAUTION: National laws and regulations must be observed.

i The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

7.1 Reprocessing of accessories

▷ The reprocessing of accessories is described in the respective manual.

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Reprocessing

7.2 Cleaning and disinfection

NOTE: Damage to the HF device due to incorrect reprocessing: No method other than the wipe-down disinfection described here is permitted.

WARNING: Fire hazard and risk from electric shock! Disconnect the plug before cleaning. Only use approved cleaning/disinfection agents for surface cleaning as per the manufacturer's instructions. It is essential to ensure that no liquids can penetrate the housing. Make sure that the AUTOSTART function is deactivated.

Wipe-down disinfection of devices and footswitches

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

KARL STORZ has qualified the resistance of the surfaces vis-a-vis wipe-down disinfection with the agent *Dr. Schumacher – Cleanisept Wipes*.

Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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8 Maintenance and repair

8.1 Maintenance

 **DANGER:** Risk of infection! Carry out a surface disinfection and wrap the device in addition to the shipping packaging material before allowing the device to leave the hospital or practice to avoid spreading germs and infections.

- ▷ Check the device, the equipment cart and the accessories (e.g. footswitch, cable) after each use for damage or defects. In particular, make sure that the insulation is intact on all cables.
- ▷ Do not use any damaged device, damaged equipment cart or damaged accessories.
- ▷ Replace defective accessories immediately.
- ▷ Have the safety inspection for the device performed once a year. Please consult and comply with the respective service manual for additional technical information.

8.1.1 Safety inspection

Safety inspections must be performed once a year.

- ▷ The next safety inspection date of the AUTOCON® III 400 can be displayed in the menu, see section Menu 'System settings', page 71.

■ A warning message appears during system start-up if a safety inspection is due. Press OK to confirm this message.

 Any shorter safety inspection cycles specified in national regulations must be observed.

- ▷ The device and accessories may be inspected only by persons who have the required training, knowledge or experience and who can perform the inspection independently.
- ▷ With regard to the safety inspection, you must comply with the country-specific rules and regulations.

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The inspector documents the inspection results and measured values corresponding to the printed inspection protocol.

In the case of severe deviations from the values of the attached final acceptance report, or if the specified maximum values were exceeded:

- ▷ Send the HF device to the service center, see section 9.1 – Technical Support, page 100.

8.2 Repair

NOTE: You can damage the HF device by doing your own repairs and modifications of medical equipment! If a repair is necessary, only allow this to be done by the service center specified (see section 9.1). Never carry out any repairs yourself.

KARL STORZ is liable for safety, reliability and performance of the HF device under the following conditions:

- Full compliance with all instructions regarding the installation and proper use for the intended purpose contained in this instruction manual was maintained.
- Changes, repairs, new settings and similar procedures were carried out only by persons authorized to do this work by KARL STORZ.
- The electrical installations in the relevant room meet the local requirements and statutory provisions.

Fast and satisfactory repairs can only be guaranteed when all required data have been supplied in full.

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The following information is required for returning the device:

- Full address
- Model number
- Serial number
- Software version

- ▷ Describe the problem, the corresponding application and the accessories used.
 - or –
- ▷ Describe the repairs to be made.

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9 Storage

- ▷ If you store the HF device for longer than one year, pay specific attention to the indicators during automatic functional testing, see section Test for proper functioning, page. 26.
- ▷ Clean the HF device thoroughly before you put it into storage.
- ▷ Store the HF device in a clean, dry place in accordance with the storage conditions.

Storage conditions:

- Temperature: -20 °C to +50 °C
- Relative humidity: 0 to 75 %, non-condensing
- Air pressure: 500 to 1060 hPa

9.1 Technical Support

In Germany contact the following service center for maintenance and repair work:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen, Germany
Service hotline: +49 7461/708 980
E-mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ subsidiary or authorized dealer.

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10 Technical data

10.1 Technical data for
AUTOCON® III 400

Insulation type / Classification	
EMC	IEC 60601-1-2
Protection class acc. to IEC 60601-1	I
Applied part type according to IEC 60601-1	CF
Standard compliance	IEC 60601-1: 2005+A1:2012, IEC 60601-1-2: 2007, IEC 60601-2-2: 2009, ISO 14971: 2007, ISO 13485: 2003 + Cor.1 2009
Classification according to Annex IX of the EC Directive 93/42/EEC	IIb
This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.	CE 0123

i The code number after the CE mark indicates the responsible notified body.

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Programs	
Number of programs in the device	300
Default programs, factory set	Yes
Individually programmable	Yes
Information shown on the display	Yes

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Neutral electrode monitoring	
EASY: Electrode Application System	Yes
Display of one-piece or split or Baby electrode	Main menu and neutral electrode selection
Contact resistance between individual sections of split neutral electrodes shown on display	Using color and contact indicator
Lead resistance shown on the display when a non-split neutral electrode is used	Yes
Maximum permissible resistance between the sections of a split electrode	300 Ω
Warning signal for hazardous conditions in connection with neutral electrodes	Visual, acoustic
Tones	Warning, activation and button tones, start melody
Warning message on the display	Text message with further information

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Safety features	
ISSys (Integrated Safety System)	Yes
Arc control	ARC CONTROL
Continuous monitoring of HF leakage current and error message	Text message with further information
Dosage monitoring with error message on the display	Yes
Continuous self-test	Yes
Continuous status indication on the display	Yes
Operating errors shown on the display	Text message with further information
System errors shown on the display	Text message with further information
Safety inspection	Automatic reminder function (optional)
Instruction manual	Hard copy

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Documentation	
Data acquisition and storage in the device	System information
Documentation of error states	Yes
Documentation of operating errors	Yes
Retrieval of system messages via the display	Text message with further information

Communication	
Display	Capacitive touch screen 9"
External interface for communication between HF generator and an Argon beamer	Fiber-optic cables
USB interface for software updates	Yes
External PC interface for service support using KARL STORZ software	UART
External interface for communication between the HF generator and the KARL STORZ OR1™ system.	Ethernet connection

Service support	
Network connection for service support	Yes
Service support by service programs integrated in the device	Yes
Service support via ISSys	Yes

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Cooling	
Convection	Yes
Temperature-controlled fan	Yes

Operating mode	
Operating mode	Intermittent 10 s/30 s (on/off)

Specifications	
Max. monopolar power	400 W (at 200 Ω)
Max. bipolar power	400 W (at 75 Ω)
Output frequency	350 kHz / 1 MHz
Unipolar sockets	2x (footswitch and fingerswitch)
Bipolar sockets	3x (footswitch and fingerswitch)
Connection for footswitch	2x
AUTOSTART	Yes
Options	BiVascularSafe (UH 401 and UH 401U)
Scope of supply (see section 3.3)	Inspection protocol

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


Conditions of operation, transport and storage	Operation	Transport and storage
Temperature	+10 °C to +40 °C	-20 °C to +60 °C
Relative humidity	30 to 75 %, non-condensing	0 to 90 %, non-condensing
Air pressure	700 to 1 060 hPa	500 to 1 060 hPa
Operating altitude (max.)	3000 m above sea level	

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








i These max. values are not necessarily created at rated load.
The HF power limitation is subject to a tolerance of $\pm 20\%$.

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


10.2 Output, voltage and current diagrams

Unipolar cutting – laparoscopy

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Laparoscopy' = 100 W

- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Laparoscopy'.
Rated load resistance = 500 Ω

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Laparoscopy' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Laparoscopy' (idle mode)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Standard' = 200 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Standard'. Rated load resistance = 200 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Standard' = 400 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Standard' (idle mode)

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

Measurement at ohmic resistances

Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Micro' = 25 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Micro'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Micro' = 50 W

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
- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Micro' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection'

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- Table of power output P [W] as a function of the setting 'Unipolar Cutting Resection'. Rated load resistance = 500 Ω

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Resection' (idle mode)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection C-Cut®' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Resection C-Cut®'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection C-Cut®' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Resection C-Cut®' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting SupraLoop'

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- Table of power output P [W] as a function of the setting 'Unipolar Cutting SupraLoop'. Rated load resistance = 100 Ω

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting SupraLoop' (idle mode)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Dry' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Dry'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Dry' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Dry' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 1' coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Gastro Loop 1'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 1' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Loop 1' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 2' coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Gastro Loop 2'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 2' cut phase

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
- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Loop 2' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 3' coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Gastro Loop 3'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 3' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Loop 3' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 1' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Gastro Knife 1'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 1' cut phase

Effect	U (Vp)
1	650
2	650
3	650
4	700
5	750

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Knife 1' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 2' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Gastro Knife 2'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 2' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Knife 2' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 3' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Gastro Knife 3'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 3' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Knife 3' (idle mode)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Argon' = 150 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Argon'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Argon' = 300 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Argon' (idle mode)

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
Unipolare Koagulation – Laparoskopie

Unipolar coagulation – Laparoscopy

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Laparoscopy' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Laparoscopy'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Laparoscopy' = 120 W

- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Laparoscopy' (idle mode) = 1800 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Moderate' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Moderate'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Moderate' = 120 W

Effect	U (Vp)
1	250
2	250
3	250

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Moderate' (idle mode)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced coag' = 40 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced coag'. Rated load resistance = 1000 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced coag' = 80 W

- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Forced coag' (idle mode) = 3500 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Resection' = 60 W

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- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Resection'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Resection' = 120 W

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Spray' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Spray'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Spray' = 120 W

Effect	U (Vp)
1	3000
2	3800
3	4600
4	5000

- Diagram of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Spray' (idle mode)

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
Unipolare Koagulation – Forciert mixed

Unipolar coagulation – Forced mixed

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced mixed' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced mixed'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced mixed' = 120 W

Effect	U (Vp)
1	1500
2	2000
3	2500

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Forced mixed' (idle mode)

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
Unipolare Koagulation – Forciert cutting

Unipolar coagulation – Forced cutting

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced cutting' = 125 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced cutting'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced cutting' = 250 W

Effect	U (Vp)
1	1500
2	1500
3	1300
4	1300


- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Forced cutting' (idle mode)

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Unipolare Koagulation – Gastro Coag

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Gastro Coag' = 25 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Gastro Coag'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Gastro Coag' = 50 W

Effect	U (Vp)
1	1800
2	2200
3	2800

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Gastro Coag' (idle mode)

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
Unipolar coagulation – Argon flexible
(only in combination with an Argon beamer)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flexible' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Argon flexible'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flexible' = 120 W

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
- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Argon flexible' (idle mode) = 4400 Vp

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
Unipolar coagulation – Argon flex. pulse
(only in combination with an Argon beamer)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flex. pulse' = 40 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Argon flex. pulse'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flex. pulse' = 80 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Argon flex. pulse' (idle mode)

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i The pulse frequency changes with the effect setting. The higher the effect level, the faster the pulse sequence.
Effect 1: 1 Hz, effect 2: 5 Hz, effect 3: 10 Hz
The mode 'Argon flexible' is paused due to these pulse sequences.

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Unipolar coagulation – Argon open
(only in combination with an Argon beamer)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon open' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Argon open'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon open' = 120 W
- HF output voltage U [Vp] for the setting 'Unipolar Coagulation Argon open' (idle mode) = 4600 Vp

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Unipolare Koagulation – Cardiac Thorax

Unipolar coagulation – Cardiac Thorax

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Thorax' = 50 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Cardiac Thorax'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Thorax' = 100 W

- HF output voltage U [Vp] for the setting 'Unipolar Coagulation Cardiac Thorax' (idle mode) = 1800 Vp

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Unipolare Koagulation – Cardiac Mammaria

Unipolar coagulation – Cardiac Mammary

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Mammary' = 30 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Cardiac Mammary'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Mammary' = 60 W
- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Cardiac Mammary' (idle mode) = 1800 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation SimCoag' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation SimCoag'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation SimCoag' = 120 W

Effect	U (Vp)
1	2000
2	2500
3	4600

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation SimCoag' (idle mode)

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Cutting' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Cutting Bip. Cutting'.
Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Cutting' = 200 W

- HF output voltage U [Vp] with the setting 'Bipolar Cutting Bip. Cutting' (idle mode) = 400 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bipolar Scissors' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Cutting Bipolar Scissors'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bipolar Scissors' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Cutting Bipolar Scissors' (idle mode) = 200 Vp

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Resection' cut phase

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
- Table of output power P [W] as a function of the setting 'Bipolar Cutting Bip. Resection'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Resection' phase after cutting

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
- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Cutting Bip. Resection' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Vaporization'

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- Table of power output P [W] as a function of the setting 'Bipolar Cutting Bip. Vaporization'.
Rated load resistance = 75 Ω

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- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Cutting Bip. Vaporization' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation RoBi®' = 50 W

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

- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation RoBi®', Rated load resistance = 25 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation RoBi®' = 100 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation RoBi®' (idle mode) = 110 Vp

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Measurement at ohmic resistances



- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Laparoscopy' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Laparoscopy'.
Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Laparoscopy' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Laparoscopy' (idle mode) = 150 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard' = 120 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard' (idle mode) = 150 Vp

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
Bipolare Koagulation – Bip. Resektion

Bipolar coagulation – Bip. resection

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bip. Resection' = 175 W

Effect	P (W)
1	125
2	200
3	275
4	350

- Table of power output P [W] as a function of the setting 'Bipolar Coagulation Bip. Resection'. Rated load resistance = 75 Ω

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- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Coagulation Bip. Resection' (idle mode)

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
Bipolare Koagulation – Bip. Vaporisation

Bipolar coagulation – Bip. vaporization

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bip. Vaporization'

Effect	P (W)
1	250
2	250
3	250

- Table of power output P [W] as a function of the setting 'Bipolar Coagulation Bip. Vaporization'.
Rated load resistance = 75 Ω

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
- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Coagulation Bip. Vaporization' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'BiVascularSafe'/'BiVascularSafe Open Surgery'

- Power output P [W] with the setting 'BiVascularSafe'/'BiVascularSafe Open Surgery' (Rated load resistance = 25 Ω) = 200 W

- HF output voltage U [Vp] with the setting 'BiVascularSafe'/'BiVascularSafe Open Surgery' (idle mode) = 190 Vp

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

Bipolar coagulation – Standard AUTO

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard AUTO' = 60 W


- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard AUTO'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard AUTO' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard AUTO' (idle mode) = 150 Vp

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Bipolare Koagulation – Bipolare Schere

Bipolar coagulation – Bipolar scissors

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bipolar Scissors' = 60 W

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


- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Bipolar Scissors'.
Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bipolar Scissors' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Bipolar Scissors' (idle mode) = 200 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Micro' = 10 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Micro'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Micro' = 20 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Micro' (idle mode) = 90 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Forced' = 50 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Forced'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Forced' = 100 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Forced' (idle mode) = 550 Vp

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10.3 Technical documentation

On request, the manufacturer will provide those service manuals, itemized parts listings, descriptions, sets of adjustment instructions and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

 We reserve the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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11 Accessories/Replacement parts

Original KARL STORZ accessories are suitable for use with the AUTOCON® III 400 devices. When using accessories made by other manufacturers, the user must ensure that they are designed for and compatible with the maximum HF peak voltage of the HF device.

For the use and correct reprocessing of the autoclavable devices, compliance with the relevant instruction manuals accompanying these devices is required.

Detailed information on accessories and replacement parts is available in the corresponding KARL STORZ catalogs.

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12 EMC

12.1 Guidance and manufacturer's declaration as per IEC 60601-1-2, section 6.8.3.201

Electromagnetic interference emissions (IEC 60601-1-2, Table 1)		
AUTOCON® III 400 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 400 must make sure that it is used in such an environment.		
Interference emissions measurements	Compliance	Electromagnetic environment – guidelines
HF emissions acc. to CISPR 11	Group 2	The AUTOCON® III 400 must emit electromagnetic energy in order to perform its intended function. Nearby electronic devices may be affected.
HF emissions acc. to CISPR 11	Class A	The AUTOCON® III 400 is solely intended for use by qualified medical staff. The AUTOCON® III 400 can cause radio interference, e.g. with radio receivers, or may interfere with the operation of equipment in the vicinity. It may be necessary to introduce appropriate corrective measures, e.g. re-orienting, relocating, or shielding the unit or filtering the connection to the public mains network.
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	

Electromagnetic immunity (IEC 60601-1-2, Table 2)

The AUTOCON® III 400 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 400 must make sure that it is used in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	±6 kV contact discharge	±6 kV contact discharge	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity must be at least 30 %.
	±8 kV air discharge	±8 kV air discharge	
Rapid transient electric interferences/ bursts acc. to IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
	±1 kV for input and output lines	±1 kV for input and output lines	
Surges acc. to IEC 61000-4-5	± 1 kV outer conductor – outer conductor	± 1 kV outer conductor – outer conductor	Mains power quality should be that of a typical commercial or hospital environment.
	± 2 kV outer conductor – earth	± 2 kV outer conductor – earth	
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	$< 5 \% U_T$ for $\frac{1}{2}$ cycle (> 95 % dip) $40 \% U_T$ for 5 cycles (60 % dip) $70 \% U_T$ for 25 cycles (30 % dip) $< 5 \% U_T$ for 5 s (> 95 % dip)	$< 5 \% U_T$ for $\frac{1}{2}$ cycle (> 95 % dip) $40 \% U_T$ for 5 cycles (60 % dip) $70 \% U_T$ for 25 cycles (30 % dip) $< 5 \% U_T$ for 5 s (> 95 % dip)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AUTOCON® III 400 requests continued functioning even when interruptions of the power supply occur, it is recommended that the AUTOCON® III 400 be powered from an uninterruptible power supply or a battery.

Note: U_T is the a. c. mains voltage prior to application of the test level.

Electromagnetic immunity – for equipment and systems that are not life-supporting (IEC 60601-1-2, Table 4)

The AUTOCON® III 400 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 400 must make sure that it is used in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Conducted HF disturbances acc. to IEC 61000-4-6 Radiated HF disturbances acc. to IEC 61000-4-3	3 V_{effective value} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V_{effective value} 3 V/m	Portable and mobile HF communications equipment should be used no closer to any part of the AUTOCON® III 400, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. The field strength of stationary transmitters at all frequencies on site ^a should be, according to a study, less than the conformity level ^b. Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.		
a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed HF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AUTOCON® III 400 is used exceeds the applicable HF compliance level above, the AUTOCON® III 400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the AUTOCON® III 400.		
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended separation distances between portable and mobile HF communications equipment and the AUTOCON® III 400 (IEC 60601-1-2, Table 6)

The AUTOCON® III 400 is intended for use in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the AUTOCON® III 400 can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF communications equipment (transmitters) and the AUTOCON® III 400 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance d according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
2	1.7	1.7	3.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.		

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13 Disposal

 The national regulations must be observed when disposing of or recycling the device or its components.

Symbol	Designation
	Devices bearing this symbol must be collected separately as electrical and electronic devices. Within the EU, disposal is handled free of charge by the manufacturer.

▷ For questions on the disposal of the device, please contact the Technical Support, see section Technical Support on page 100.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kse.com

KARL STORZ Veterinary Endoscopy-America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098
E-Mail: karlstorz@enet.cu

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street

2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building West Side 7A - Unit 7WA - 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
14th Floor, Kumho Asana Plaza
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 (8)3832 8000, Fax: +84 (8)3832 8039
E-Mail: karlstorz@hcm.vnn.vn

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Onon Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Janguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: stephen.Lee@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeonggil-Ho, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw#

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latin America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopy
Miramar Trade Center
Edificio Jerusalem. Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapulterc Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America do Sul Ltda
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, São Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsatragrand 14
12724 Skarholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4. Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopie Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopie France S.A.S.
12, rue Georges Guynemer. Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: stolz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) - Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-I
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 - Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: ningwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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Important information for users of KARL STORZ instruments

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the AUTOCON® II 400. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

NOTE: The instruction manual only applies to the following unit versions:
20 5352 20-122/20 5352 20U122
20 5352 20-125/20 5352 20U125
20 5352 21-125/20 5352 21U125
20 5352 22-125/20 5352 22U125
20 5352 23-125/20 5352 23U125

(See Section 'Unit versions', p. 62).

General warning

Users of AUTOCON® II 400 devices should be thoroughly trained in its use and the applicable medical procedures. Use of this equipment and accessories should not be undertaken until the user has fully familiarized himself/herself with the instructions for use, assembly and care. Instruction manuals should be carefully studied and be available to the surgical team during the procedure; it is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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Controls, displays, connectors, and their uses

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- ① Power switch, illuminated
0 = off
- ② Touch screen
- ③ HF output sockets/neutral elec.; socket
assignment depending on unit versions
(see p. 62)
- ④ Neutral electrode
- ⑤ Connection socket for footswitch
(3-pedal footswitch 20017831 or
2-pedal footswitch 20017830 and
1-pedal footswitch 20017832)
- ⑥ Service connection for software updates
- ⑦ Connectors for SCB (for **SCB**® see page VI)
- ⑧ Potential equalization connector
- ⑨ Non-heating appliance coupling socket
with line fuses 2 x T4.0 AH250V at 230 V or
2 x T8.0 AH250V at 115 V
- ⑩ Resection auxiliary module 20535280 (con-
nect the auxiliary module to the multifunc-
tion socket on the AUTOCON® II 400 model
2053522x-(U)125 using the integrated cable.
Connect the instrument to the Bi-Saline resec-
tion socket on the auxiliary module using the
bipolar cable 27176 LEB/LEBL; see photo
on page IV.) (Please also refer to the instruc-
tion manual 96206181 DER for the auxiliary
module.)
- ① **NOTE:** The auxiliary module 20535280
is included in the scope of supply for the
AUTOCON® II 400 models 20535221-(U)125
and 20535223-(U)125.

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Service connection for software updates

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Potential equalization connector

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Applied part type CF as per IEC 60601-1 (the
neutral electrode may remain applied to the patient
during defibrillation)

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The SCB interface, based on the CAN field bus,
permits remote control of equipment functions, as
well as remote display of equipment parameters.

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Environmental protection use period of 50 years
(China RoHS)

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Dual/triple-pedal footswitch approved for cutting only (yellow footswitch)

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Dual/triple-pedal footswitch approved for coagulation only (blue footswitch)

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Yellow footswitch button for cutting and blue footswitch button for coagulation

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Single-pedal footswitch – only for coagulation

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+/- button for setting the cutting and coagulation parameters: power and effect

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Indicator lights green when button is selected

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General warning	III
Images of the equipment	IV
Controls, displays, connectors and their uses	V
Symbols employed	VI

General information	1
The principle of high frequency surgery	1
The AUTOCON® II 400 unit	7

Safety instructions	8
Warnings and cautions	8
Intended use	13
User qualification	13
Patient profile	13
Safety precautions at the site of installation	15
Safety precautions when operating the system	15
Important notes on the use of the unipolar technique	16
Safety devices	23

Definitions and explanations for AUTOCON® II 400	31
Definitions and explanations	31
Symbols for effects which can be set	59
Unit versions	62

Installation and operating instructions	64
Unpacking the equipment	64
Basic equipment	64
Installing and connecting up the unit	64
Operating the unit	67

Maintenance	89
Fuse replacement	89
Reprocessing	90
Maintenance and safety check	92
Maintenance	92
Safety check	92
Servicing and repair	93
Disposal	93
Repair program	94
Limitation of liability	95
Manufacturer's warranty	95

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**Spare parts,
recommended accessories..... 130**

Spare parts list.....	130
Recommended accessories.....	130

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The principle of high frequency surgery

Biological tissue contains a more or less high concentration of electrolytes, making it sufficiently conductive to be treated electrosurgically. The thermal effect of HF current is used for separating (cutting) and coagulating tissue (desiccation of tissue and/or hemostasis).

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
HF currents must be used on the patient since low-frequency currents can stimulate nerve and muscle cells in the current flow due to electrochemical processes (electrolysis). These effects are small enough to be disregarded with frequencies above 100 kHz. In HF surgery, a distinction is made between the following application techniques:

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
Unipolar operation techniques

With the unipolar operation technique, the electric current flows from the small-face active electrode to a large-face passive or neutral electrode, the patient's body always forming part of a closed circuit. Tissue separation or coagulation takes place at the active electrode.

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
Bipolar operation techniques

With the bipolar operation technique, the current flows between two paired electrodes, the faces of which are either the same or different sizes, depending on the application (cutting or coagulation). Only a small part of the tissue is included in the circuit between the electrodes during application.

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Separation of tissue

When separating tissue (cutting), the cell sap of the tissue is vaporized. Due to the resulting vapor pressure, the tissue cells tear and the escaping vapor pressure causes the electrode to be insulated. If the voltage is high enough, the current flow is maintained by the commencement of arcing between the electrode and tissue. The thermal effect of the current causes the separation of tissue and the formation of surface slough, depending on the size of the arc.

Arcing commences at voltages of around 200 V~. If voltages above 500 V~ are used, the result is substantial formation of slough due to the energy released in the arc.

For cutting, automatic control loops allow an optimal, reproducible cutting quality to be achieved, which is to the greatest possible extent independent of cutting depth, cutting speed, tissue resistance and such like. Power output is automatically adjusted to the level required for cutting.

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Tissue coagulation

To coagulate biological tissue, it must be heated to approx. 70 °C.

To prevent the current from causing tissue separation, the power supply, and consequently the amount of heat, must be used specifically to produce denaturation of the tissue.

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The thermal effects which take place in biological tissue, depending on temperature, can be described as follows:

- Up to approx. 40 °C: No visible cell damage
- Approx. 40 °C and above: Reversible cell damage depending on the duration of current flow
- Approx. 49 °C and above: Irreversible cell damage
- Approx. 70 °C and above: Conversion of collagens into glucose – shrinking of tissue containing collagen – hemostasis
- Approx. 100 °C and above: Dehydration/desiccation: phase transition of intracellular and extracellular water from liquid into gas. Adhesive effect of glucose after dehydration. Shrinking of the coagulum after dehydration.
- Approx. 200 °C and above: Carbonization and fourth-and above: degree pathological burns. Unpleasant odor of burnt tissue.

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Due to the inhomogeneity of the electrical and thermal properties of tissue and the resulting uneven current density distribution within the tissue, the temperature rises with an irregular distribution when power is supplied.

The current density, and thus the temperature, is highest in the vicinity of the coagulation electrode. The temperature decreases very rapidly with increasing distance from the electrode so that coagulation is at first limited to the area surrounding the coagulation electrode.

The effective contact area between the tissue and coagulation electrode has a decisive influence on the maximum achievable area of the coagulation zone.

Due to the spatial extent of the coagulation zone being dependent on the current density distribution, the technique with which the HF current is applied is also of importance. Here, a distinction is made between unipolar and bipolar surface application, unipolar and bipolar puncture application as well as unipolar and bipolar grip application.

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- Standard-Koagulation
- Forcierte Koagulation
- Spray-Koagulation

Coagulation modes

The spatial extent of the coagulation zones depends on the type of energy supply and application technique.

It is useful to differentiate between 3 coagulation modes:


- Standard coagulation
- Forced coagulation
- Spray coagulation

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
Standard coagulation

With standard coagulation, the generator voltage is kept so low that no arcing is able to take place between the electrode and tissue during the coagulation process. Standard coagulation is recommended for all instances of coagulation in which unipolar or bipolar coagulation electrodes are used which are brought into direct contact with the tissue that is to be coagulated. Standard coagulation provides for deep coagulation while avoiding carbonization of the tissue.

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Forced coagulation


With forced coagulation, electric arcs are intentionally permitted between the coagulation electrode and tissue to enable the formation of surface sloughing with nevertheless deep-reaching coagulation. This particularly applies if only small-face electrodes are available for coagulation. To avoid the risk of a cutting effect, a modulated HF voltage must be used. Tissue sloughing will, however, occur.

A typical application is TUR, in which the cutting electrodes are also used for coagulation purposes.

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
Spray coagulation

If sufficiently high voltage pulses are used for coagulation, arcs occur between the electrode and tissue of such a length that contact between the tissue and the electrode is no longer necessary. Spray coagulation is primarily used for surface coagulation or for hemostasis of vessels which cannot be directly reached with a coagulation electrode (e.g. in bone fissures).

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Legend:

- A: Power required for cutting
B: Power output without control of unit
C: Excess power

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The operator bases his assessment of the quality of the tissue cutting on the various parameters of the entire system.

Such parameters include the geometry and size of the electrode used, the nature of the tissue at the site of the operation, the penetration depth and cutting speed of the electrode, the conductivity of the irrigant used, and the application of the neutral electrode etc.

These 'external' parameters, which can be extremely variable, influence the properties of conventional HF equipment to such a degree that the operator must constantly readjust the power or generally set a higher power level. The excess power, however, produces an unnecessary degree of surface slough formation and exposes the patient to unnecessarily high current and the associated risks.

Fig.: Power output (P) in relation to time (t) when cutting tissue.

The cutting action during tissue cutting is influenced to a great extent by the electric arc maintained between the electrode and the tissue. This realization led to the development of a new generation of high frequency surgical units.

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The AUTOCON® II 400 unit

The AUTOCON® II 400 HF surgical unit is a universally applicable unit which permits the precise setting of various cutting and coagulation qualities as well as reliable reproducibility of these settings. Due to automatic closed-loop control of the HF voltage, the HF voltage and intensity required for achieving the desired cutting quality are provided in each phase of the incision.

The AUTOCON® II 400 permits both unipolar and bipolar incisions and provides a suitable cutting and coagulation mode for every application.

The bipolar coagulation mode additionally permits operation with the autostart and autostop functions.

In conjunction with automatic closed-loop control of cutting and coagulation quality, the various safety devices of the instrument provide an exceptionally high level of user and patient safety. Software-supported test programs ensure easy and rapid servicing.

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Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be read carefully to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

Definitions

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

NOTE: A Note indicates special information about operating the unit.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

WARNING: Keep out of reach of patients.

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

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Safety instructions

Warnings and cautions

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

WARNING: To avoid the risk of an electric shock, this unit may only be connected to a power supply network with a protective conductor.

WARNING: Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly. The unit should not be used if any damage is evident.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: Always unplug the unit before carrying out any maintenance work on it.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Contact with noninsulated instruments by the active high frequency electrode can lead to unintentional burns. We therefore recommend using insulated instruments in the vicinity of the active electrode.

WARNING: Avoid currents flowing via tissue bridges. The resulting current concentrations at the base adjoining healthy tissue may cause burns.

WARNING: Use standard-compliant endoscopes with an insulated eyepiece to reduce the risk of capacitive leakage currents flowing to the user or to the attached camera system.

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Safety instructions

Warnings and cautions

WARNING: Subjecting the instrument to too high electrical loads can damage it. Should the damaged area come into contact with tissue, it may result in undesired coagulation.

Ascertain the electrical load capacity of the instrument. This will either be printed on the instrument or found in the corresponding instruction manual. Compare the electrical load capacity of the instrument with the maximum HF peak voltage of the desired mode (see section "Technical data").

WARNING: The arc produced at the electrode during operation of the high frequency surgical device can impair the function of other electrical devices.

WARNING: During simultaneous use of high frequency surgery and physiological monitoring equipment on a patient, every monitoring electrode should be positioned as far away as possible from the active and neutral electrodes of the high frequency surgical device. Needle electrodes for monitoring are not recommended.

In all cases monitoring electrodes which contain devices for limiting the high frequency current are recommended.

WARNING: A double-faced neutral electrode that has short-circuited precludes monitoring.

WARNING: The unit may only be operated under the environmental conditions specified.

WARNING: Please ensure that there is sufficient air circulation for your unit.

WARNING: Please always observe cleaning and sterilization instructions.

WARNING: The unit must not be operated if the touch screen or display is defective

WARNING: It is not permitted to press more than one touch point simultaneously.

WARNING: Patient contact with grounded metal parts must be avoided.

WARNING: In the event of a sudden power loss, check that the neutral electrode is correctly positioned before increasing the power output.

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Safety instructions

Warnings and cautions

WARNING: Due to low-frequency electric currents, nerves and muscles of the patient can be stimulated during HF surgery.

WARNING: Patients with heart pacemakers must be treated with special care. High frequency energy can interfere with heart pacemakers. The pacemaker manufacturer's instructions must be strictly followed.

WARNING: Failure of the HF surgical equipment could result in an unintended increase of output power.

Use of several HF units on one patient

WARNING: Only use units with 'floating' application parts. Always operate a unit with its own neutral electrode.

Avoid high frequency electric circuits, i.e. ensure adequate spacing (i.e. at the top, bottom, in the center, no switching of active to neutral electrodes). All units should be operated according to the instructions provided in the respective instruction manuals.

Injury due to incorrect high frequency parameter

WARNING: The mode and power setting should be adapted to the intervention and should be as low as possible. This avoids the following:

- thermal invasion of the tissue, which can be caused by too high a power or effect intensity;
- insufficient coagulation due to too low a power setting, which can result in excessive bleeding.

WARNING: When performing unipolar HF surgery on bio-tissue with fine structures, there is the risk of burns due to local current concentration. For this reason, bipolar HF surgery is recommended in these areas.

WARNING: Avoid explosive gas concentrations that may be present in the area where a high frequency surgical procedure is being conducted. The following situations, among others, must be taken into consideration:

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Safety instructions

Warnings and cautions

Gas which may support combustion is sometimes present in the gastrointestinal tract of an unprepared patient. Certain substances used for the preparation of the lower GI tract can enhance methane production. This is particularly relevant in colonoscopy, but has also been recorded in the upper GI tract. In addition, during transurethral resection of the prostate, it has been recorded that hydrogen can accumulate in the bladder above the irrigation solution, this is particularly true when a strong arc is used. Working with explosive gas mixtures may lead to deflagration/explosion.

WARNING: Use only fuses of the correct rating.

WARNING: Avoid allowing liquids to enter the unit. Do not store liquids on or above the unit. It must also be ensured that no liquid bottles or tubing are placed near the unit.

If liquid enters the unit, despite the precautions, sufficient time should be planned for evaporation (this also applies for the formation of condensation, see "acclimatization", p. 120).

CAUTION: Only operate the unit with the voltage stated on the unit's identification plate.

CAUTION: The signal input and signal output parts of this equipment are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

CAUTION: It is essential to follow the instruction manual closely when cleaning the unit and preparing/sterilizing unit accessories.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

NOTE: In the case of HF generator accessories, attention must be paid to the maximum HF peak voltage allowed. Please observe the instruction leaflet or instruction manuals provided to prevent damage to the material.

NOTE: Inspect high frequency cord, the connecting cable to the neutral electrode and endoscopic accessories regularly for damage to the insulation.

NOTE: KARL STORZ recommends using a smoke evacuation system for the removal of smoke generated during the use of an HF surgical unit.

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Before using the unit on the patient for the first time it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

High frequency surgical devices are universally applicable treatment units for unipolar and/or bipolar cutting and/or coagulation of biological tissue.

The AUTOCON® II 400 may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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User qualification

The AUTOCON® II 400 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the instrument. Such instruction should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and information on further training alternatives.

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Patient profile

The use of AUTOCON® II 400 is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

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Affected body part or tissue type

The field of application is not confined to certain parts of the body or types of tissue.

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User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Adequate knowledge of the purpose and whereabouts of the medical products logbook
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

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Use

The device may only be used in hospitals and doctors' offices in the ambient conditions given in the technical data (page 120).

Other conditions of use:


- Frequency of use up to several times a day
- Length of use from a few minutes to several hours a day
- Position it on a vibration-free, fixed base
- Can be moved only by installing it on an equipment trolley
- Can be swiveled only by installing it on ceiling supply units
- Intermittent operation:
25% duty cycle (10 sec. on, 30 sec. off)

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Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones.

This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected in accordance with national regulations.

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Safety precautions when operating the system

It is the user's responsibility to make sure the unit is safe and operates properly before using it.

During treatment with the AUTOCON® II 400 the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if visual observation of the action of the unit is ensured.

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Important notes on the use of the unipolar technique

Risk of burns

HF surgery always entails certain risks for both user and patient. Above all, contact burns can occur due to unintentional thermal damage of tissue. With a knowledge of the possible causes, and by observing the safety instructions below, burns can, however, be avoided.

Incompatible or single-faced neutral electrodes

When using an incompatible neutral electrode, you must take into account that the monitoring of the contact between the neutral electrode and the skin is not accurate. When using a single-faced neutral electrode, the contact between the neutral electrode and the skin is not monitored. The device does not produce an optical or acoustic warning signal when there is insufficient contact between the neutral electrode and the skin. This exposes the patient to a risk of burning in the area beneath the neutral electrode!

Use the accompanying documentation supplied by the neutral electrode's manufacturer to establish whether the neutral electrode is suitable for the AUTOCON® II 400 in use.

Only suitable neutral electrodes should be used.

When using a single-faced neutral electrode: Check regularly that the neutral electrode is in appropriate contact with the skin.

Use the accompanying documentation supplied by the neutral electrode cable's manufacturer to establish whether the neutral electrode cable is suitable for the neutral electrode in use.

Only suitable neutral electrode cables should be used.

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Burns due to incorrect application of the neutral electrode

With incorrect or no application of the neutral electrode, there is a high risk of thermal tissue damage both at the point of application of the neutral electrode and at other points of the patient's body. When applying neutral electrodes, it is therefore imperative that care be taken to ensure that the effective contact area or conductance between the neutral electrode and the patient is adequate for the HF power used or the intensity of the HF current. Here, effective contact area refers to the area of the neutral electrode in electroconductive contact with the patient's skin during HF surgery.

When selecting the neutral electrode, take the following empirical values into account:

Neutral electrode made of electroconductive silicone:	1.5 watts/cm ² or 10 mA _{rms} /cm ²
Neutral electrode made of metal without conductive jelly:	3.0 watts/cm ² or 15 mA _{rms} /cm ²
Disposable neutral electrode:	5.0 watts/cm ² or 20 mA _{rms} /cm ²

! WARNING: Care must be taken to position the neutral electrode near the operating field on shaven tissue with good blood circulation.

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
The neutral electrode is secured to either the upper arm or thigh, depending on which is nearest the operating field. To ensure that the entire face of the neutral electrode is applied, it must be securely fastened with rubber straps and then wrapped with a bandage so that it cannot slip, even if the patient should move (note the risk of decubital necrosis). The neutral electrode must always be in perfect condition. Grease, oxidation, dirt etc. on the contact surface reduces the size of the electroconductive area and can lead to dangerously high current densities.

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! WARNING: A short circuit fault is detected when using a double-faced NE. Power is disconnected until optimal conditions return. Subsequently, a renewed activation is essential.

! WARNING: The attachment of a neutral electrode made of conductive silicone (27805) prevents monitoring.

i NOTE: The HF current flows to the neutral electrode through the tissue. It should be remembered that the distribution of current at the proximal corners and edges to which the current flows may be higher than at the distal corners and edges. For this reason, when arranging the neutral electrode care must be taken to ensure that the longer side of the neutral electrode is facing the operating site.

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
Patients must be correctly positioned before performing any operations using HF surgical instruments. An electrically insulating, thick and dry pad must be placed between the patient, operating table and clamps. If moisture, perspiration etc. are anticipated during the operation, penetration of moisture through this pad must be prevented by waterproof sheeting or other suitable measures. The patient's body tissue must not be allowed to come into contact with nearby electrically conductive objects (examples: operating table, metal clamps, wet towels, etc.).

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Since contact areas will generally be small in such cases, burn-necroses could arise from high frequency leakage currents due to high local current densities. In cases where contact between the neutral electrode and patient is poor this effect is increased.

Areas subject to more severe perspiration, extremities resting against the body or skin-to-skin contact should be dried by placing cover cloths in-between. Urine should be drained off by a catheter. If an injection cannula is used as an ECG electrode, the metal cone must not rest on the skin; the same applies to the leads of monitoring equipment.

A Electrically insulated pad
B Grounded operating table

WARNING: Skin-to-skin contact should be avoided as this may lead to undesired diversions of the HF current. This incurs a risk of burning.

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
Confusing the contacts on the unipolar sockets.

CAUTION: Confusing the sockets will damage the device.

If you are using a connecting cable with a Ø 4 mm unipolar plug, the plug must only be connected to the socket with the blue ring. The correct socket is identified by an arrow in the diagram.

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Burns due to unintentional activation of the HF generator

Patients can also suffer burns due to unintentional activation of the HF generator. Unintentional activation can occur, for example, if the footswitch is inadvertently depressed, if there is a short circuit in one of the cables, if electrically conductive liquids enter the electrode handle or the footswitch, or if the HF generator is defective.

- WARNING:** To prevent such occurrences:
- never leave electrodes lying on or near the patient in such a way that could result in direct or indirect (through the mediation of electrically conductive objects) contact between electrode and patient,
 - ensure that the acoustic warning signal indicating the operational status of the HF generator is always set loud enough to be easily heard.
 - ensure that the visual warning signal can be seen by all attending personnel.

Burns from contact with hot electrodes

Cutting and coagulation electrodes can become very hot during cutting or coagulation procedures. Burns can result from active or inactive electrodes coming into contact with tissue immediately following cutting or coagulation procedures.

- WARNING:** Electrodes which are not currently in use must not come into contact with the patient.

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Burns due to improper setting of HF power output

The optimal HF power level setting is that which yields the best results in terms of the desired effect (cutting/coagulation) with minimal incidence of undesirable side effects (i.e. carbonization of tissue or sparking). Since undesirable side effects generally occur at excessively high power levels and/or if effects are too high, HF power levels should always be kept as low as possible.

If the desired effects should suddenly no longer be attainable at power levels which experience has shown to be optimal, always check to be certain that everything is in proper working order before increasing the power level:

- Is the neutral electrode properly connected?
- Is the active electrode properly connected?
- Are active and neutral electrodes clean?
- Is the HF generator power level and/or effect correctly set?

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Burns arising from ignition of combustible media

Skin cleansing agents, disinfectants, and anesthetic gases, as well as endogenic gases, may be readily combustible and can even explode. Since an arc usually occurs in HF surgery, at the active electrode, vapors or gases evolved by such media could be ignited and result in burn injuries to patients.

Thus, before switching on the HF generator, make certain that all combustible media and all objects that have been in contact with them (tissues, swabs, etc.) have been removed from the operating area.

WARNING: During TUR, when cutting or coagulating in irrigation liquid, H_2O molecules may be split into H_2 and O_2 by the electric arc. These gases may collect at the roof of the bladder in the form of a highly explosive gas mixture. Resection in this gas mixture can lead to dangerous explosions.

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Burns due to unsuitable application

If, during unipolar application, HF current has to flow through tissue constrictions, it must be remembered that a coagulation effect can occur at the constricted point. If tissue sections have a constricted cross section over an extended distance, a uniform coagulation effect can occur over the entire distance.

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Stimulation of nerves and muscles

Due to low-frequency electric currents, nerves and muscles of the patient can be stimulated during HF surgery. The arcs between the electrode and tissue, which are unavoidable for cutting processes, spray coagulation and standard coagulation, have the effect of rectifying part of the HF alternating current, whereby low-frequency current components occur which are able to cause nerve stimulation. When using HF surgery on electrically stimutable structures, contractions of the muscles concerned must therefore be expected.

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Heart pacemakers

Irreparable damage to heart pacemakers and interference with pacemaker operation, which could lead to ventricular fibrillation, are likely to result if high frequency surgery methods are used on patients with implanted heart pacemakers or pacemaker electrodes.

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Implants

In the case of patients with electrically conductive metal implants it must be ensured that with all applications involving HF surgery the current flow or current path of the HF electric current does not pass through or over the implant. At the current entry and exit points into and from the bio-tissue to the implant, current concentration may cause irreparable internal burns.

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Safety devices

WARNING: Safety circuit switches off in case of output error caused by the unit

Deviations from the setpoint of the respective HF output that has been set are corrected by the unit. An independent safety circuit monitors the power output and switches it off if proper cutting or coagulation is no longer ensured.

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Monitoring the activation time

With patient application, the HF generator is only activated briefly when used as prescribed. The activation time of the HF generator is therefore monitored. If the maximum activation time is exceeded, this is indicated by the error message A-1C on the display. The power output is then switched off. This safety circuit reduces damage resulting from unintended activation of the electrodes. The maximum activation time is limited to 99 seconds on delivery. Under 'Application Settings' you can select 'unlimited' or a specific value between 0 and 99 seconds.

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NOTE: Continuous activation can cause irreparable damage to the unit.

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Safety Devices

Monitoring the neutral electrode

The AUTOCON® II 400 has a safety circuit which monitors the electrical connection between the neutral electrode and unit and, when using double-faced electrodes, the correct application on the patient. The neutral electrode (NE) status display indicates how well the NE has been positioned on the patient. If it is positioned correctly on the patient a green neutral electrode is displayed. If a critical limit value is too low or too high, this is displayed in red.

As a general rule, the quality of the neutral electrode is indicated by its color. This may cause confusion in purely bipolar interventions (red = caution). With the exception of the menu for NE detailed display, the red NE symbol is thus suppressed when the following prerequisites are met:

- a bipolar instrument/cable (with or without resection model) is inserted in the MF socket and
- both pedals of the dual-pedal footswitch are assigned to the MF socket.

However, as soon as the unipolar socket is activated (focus changes to this socket), this suppression is relieved. It is irrelevant whether this activation leads to the application of HF energy or is canceled by the unit for safety reasons.

WARNING: If no neutral electrode is connected, the neutral electrode symbol turns red. The restriction mentioned above applies. If this is not the case, the unit must no longer be operated in unipolar mode.

If a single-faced neutral electrode is used, the status display lights up green when the unit and the neutral electrode are correctly connected. With a single-faced neutral electrode the contact quality is monitored (NE resistance), but not the current density or asymmetry. Partial detachment of a single-faced neutral electrode always increases the current concentration under the neutral electrode and results in thermal damage to tissue due to the limited contact area.

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Safety instructions

Safety Devices

When using double-faced neutral electrodes, the electrical resistance, the current density and the symmetry between the two contact areas of the neutral electrode are measured. If the resistance is outside of a tolerance range, this is recognized by the unit and a signal is given.

If the connection is interrupted or the neutral electrode positioned poorly (NE symbol lights up red) the power output in the unipolar operating modes is blocked.

If an attempt is made to activate the HF generator in a unipolar operating mode, an acoustic alarm signal is given. At the same time, an alarm message is displayed on the screen.

Neutral electrode setup setting: **any**:
Both single-faced and double-faced neutral electrodes can be connected up. In the case of the single-faced electrodes, there are no monitoring measures as described above.

Neutral electrode setup setting: **dynamic**:
is only available for double-faced electrodes. A tolerance range is placed above the currently measured impedance, which changes dynamically – always at a lower resistance. The max. permissible impedance range of up to 120 ohms is thereby reduced. Increasing the impedance deactivates the unipolar mode until a valid impedance range is reestablished.

Increased safety for patients with a low skin resistance, i.e. patients with little subcutaneous fat, children and toddlers. Even with these patients, critical detachment of the neutral electrode from the skin is detected in time.

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Safety instructions

Safety Devices

Neonatal neutral electrode monitoring

The neonatal NE monitoring is enabled when the preset neonatal procedure or a procedure created on the basis of that procedure is used. If a current limit of 300 mA is exceeded, an alarm message is shown on the display:

'A-B2: NE current limit exceeded (>300 mA)'

If the current limit is exceeded, this can be a sign that the neutral electrode is overheating. Check the neutral electrode for overheating and reduce the effect and power settings as necessary.

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
	Unit-NE connection	Skin-NE contact	NE direction of application	Added safety for patients with a low skin resistance
Double-faced NE/ dynamic NE setting	•	•	•	•
Double-faced NE/ double-faced NE setting	•	•	•	
Double-faced NE/ any NE setting	•	not entirely, pay attention to warnings	not entirely, pay attention to warnings	
Single-faced NE/ any NE setting	•			
Single-faced NE/ single-faced NE setting	•			

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Safety Devices

Self-test program

Each time the unit is powered up, a self-test program is run. First the logo is displayed for approx. 5 seconds, afterwards a white screen appears for approx. 2 seconds. This serves to optically assess if the screen is functioning correctly. If any part of the white page shown displays another color, the screen is defective and the unit should no longer be used. This is for safety reasons as the HF output data may not be displayed or displayed incorrectly. Please notify your KARL STORZ service partner at once!

If the unit detects a fault, a visual and, depending on the classification, an acoustic warning is given. The codes can be found in the table 'AUTOCON® II 400 Alarm and information messages' on page 97.

After the self-test, a procedure overview appears on the monitor. The date and time of day are displayed on the bottom right. The flashing seconds display on-screen shows that the function is activated (alive signal).

Switch-off option for HF activation

HF activation can be switched off for every HF socket and every mode by selecting 'Cut off' or 'Coag off'. Activation using the foot pedal and/or fingerswitch is then no longer possible.

This is a safety function which prevents a cutting mode from being accidentally activated for example during surgery where only coagulation is required.

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Safety Devices

Alarm system

The AUTOCON® II 400 is equipped with an alarm system. The alarm messages are displayed on the display in color and backed up with an acoustic signal.

The visual and acoustic messages used are broken down into three priority categories:

Alarm category	Color	Flashing frequency
High priority	Red	2 Hz/0.5 sec
Medium priority	Yellow-orange	0.5 Hz/2 sec
Low priority	Turquoise	Constantly on

Depending on the priority the alarm message is backed up with an acoustic signal:

High priority: An alarm sequence consists of a series of beeps of different pitch, beginning with three beeps, each for 100 ms ON and OFF, followed by a 300 ms pause, then another two beeps for each 100 ms ON/OFF. The entire sequence is repeated after a 2 s pause, following by a 5 s pause. Then the sequence starts from the beginning again (see ①).

Medium priority: An alarm sequence consists of a series of beeps of different pitch, beginning with three beeps, each for 200 ms ON and OFF, followed by a 10 s pause. Then the sequence starts from the beginning again (see ②).

Low priority: No acoustic backup signal

① NOTE: Any low priority signal will be overwritten by an alarm signal with a higher priority.

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
The alarm system is capable of identifying the alarm with the highest priority in the status bar. In the event of more than one alarm with the same priority, only the alarm which occurred first is shown in the status bar. The alarm system classifies the alarms into those which do not permit further use of the unit and those which permit further usage.

In addition to alerts, there are also purely informative messages. These messages always appear on the user interface for just 3 seconds.

The list of all error messages relevant for the user and corresponding priorities can be found in the table 'AUTOCON® II 400 Alarm and information messages' as of page 97. Alarm messages which can only be rectified by technical service personnel are listed in the service manual.

ⓘ NOTE: The user's position must always be selected so that any alarm message can be immediately identified. The distance from the unit should not be more than 4 m.

ⓘ NOTE: The sound pressure of the audible alarms is at least 65 dBA at a distance of 1 m from the front of the unit in the horizontal direction.

ⓘ NOTE: The pitch of activation and alarm signals is in the range from 1 kHz to 3 kHz.

ⓘ NOTE: The function 'Autostop reached' is confirmed both visually and acoustically.

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**Definitions and explanations for
AUTOCON® II 400****Procedures**

Procedures are specific electrosurgical interventions. A procedure is defined by a name and designation and allocated to a memory location (Example: Dr. Schmidt – Prostate resection – Pr04). In addition to user-defined procedures, each unit also includes groups of medical specialties as predefined procedures. These start at number 100 and contain the same data as the procedures defined by the user. The number and content of the predefined procedures depend on the unit version. On delivery, only the predefined procedures are available. These cannot be changed or deleted. However, they can be hidden/shown as a group of a medical specialty.

A maximum of 100 user-defined procedures can be saved. For every procedure number, a so-called level 'B' can also be generated. Creation of this second level makes it possible to switch between the two levels of a procedure with the help of the third button of a footswitch or handpiece. Level 'B' likewise contains a full data record, as described above, which can be selected normally.

A procedure can be identified as a neonatal procedure, whereby certain additional safety measures are applied. It is only ever possible to generate one user-defined neonatal procedure from a neonatal procedure; it is always only possible to generate one user-defined non-neonatal procedure from a non-neonatal procedure. The neonatal identification also applies to any level B present. The identification always applies to the complete procedure, not individual levels.

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Definitions and explanations for AUTOCON® II 400

If a procedure B is created the first one automatically becomes level A. A level B cannot exist without a level A. If a level A is deleted, level B automatically becomes level A.

A procedure contains the following data:

- procedure number (0-99)
- name (20 alphanumeric characters)
- designation (20 alphanumeric characters)
- level identification (A or B)
- current identification, i.e. currently identified as being in use by a green lit indicator
- Neonatal identification
- identification socket designation 1, 2 and 3
- cutting mode socket 1, 2 and 3
- power limitation cutting socket 1, 2 and 3
- cutting effect socket 1, 2 and 3
- coagulation mode socket 1, 2 and 3
- power limitation coagulation socket 1, 2 and 3
- coagulation effect socket 1, 2 and 3
- activation type socket 1, 2 and 3
- Identification socket 4 (NE)

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Sockets

The unit can be equipped with four different HF output sockets. At the factory, several unit versions are provided with unipolar, bipolar, multifunction and neutral electrode sockets. Details on this can be found here in the table under item 'Unit versions' (p. 62). With the exception of the neutral electrode socket, information about the general monopolar, bipolar, or multifunction types and the equipment location is available. The user can assign a socket with the following parameters:

- cutting mode
- coagulation mode
- cutting and coagulation activation type

The MF socket adopts an exceptional position as data are transferred via this socket from the connected instrument/cable:

1. Insertion is signaled by a corresponding message, followed by a second message with the approximate number of completed applications (sterilization cycles) in gradients of 5.
2. The identification via the type (unipolar, bipolar, bipolar with resection module, Bi-Vascular-Safe) is used to display the function currently assigned to the socket on the one hand, and for restriction of the mode selection on the other hand.
3. Removal is signaled by a corresponding message.

The 'bipolar resection' auxiliary module can also be connected to AUTOCON® II 400 'High End' models. The Saline++ modes can be selected with this module as of generator control software version 10.00. Depending on the model, the hardware and/or software may need to be updated in the factory.

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Neutral electrode

- dynamic
- any
- double-faced
- single-faced

Socket for neutral electrode

The socket is used for the connection of a neutral electrode in unipolar mode.

Optionally, one of the following connectors can be connected:

- ERBE neutral electrode connector with dia. 6.35 mm;
- Neutral electrode connector with 2 pins

! NOTE: The socket is equipped with a slide, which allows the connection of the connector with dia. 6.35 mm or the connector with 2 pins depending on its position (see figure).

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Modes

A mode is characterized by different generator settings such as frequency, peak voltage, modulation, etc. A total of 11 cutting and 10 coagulation modes are available. However, each specific socket is only available for a part of these modes. The user can also assign the following parameters to a mode:

- power limitation cutting,
- cutting effect,
- power limitation coagulation
- coagulation effect.

A new mode can be selected within a procedure. A list of the available modes can be found in the table on page 35. The 'Saline' modes listed are only available for the High End models.

The 'Cut off' and 'Coag off' modes represent a special case here. Selection of these modes blocks activation via the footswitch or fingerswitch. A new mode starts with the factory settings given in the table on page 40.

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**Definitions and explanations
for AUTOCON® II 400**

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*: optional

Socket	Cutting modes	Coag modes
Bipolar	Bipolar Cut Cut off	Bipolar Soft Coag Bipolar Soft Auto Stop Coag off
Multifunction	Saline-C-Cut Saline-Time-C-Cut Power Cut Top Cut LAP-C-Cut C-Cut Cut off	Bi-Vascular-Safe* Saline-Coag Saline-Time-Coag Standard Coag Forced Coag Spray Coag Coag off
Bipolar Multifunction with resection auxiliary module	Saline-C-Cut++ Saline-Time- C-Cut++ Cut off	Saline-Coag++ Saline-Time- Coag++ Coag off
Unipolar	Power Cut Top Cut Lap-C-Cut C-Cut Gastro Cut Papillo Cut Cut off	Standard Coag Forced Coag Spray Coag Coag off
NE	Resistance measurement	

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Power limitation

The power limitation can be set up to 300 Watt. However, this maximum limit is not available in every mode.

The Autostart activation type is generally only permitted up to 50 Watt; higher settings are limited to 50 Watt, including an alarm message.

It is not possible to specify a power limit in Saline and Bi-Vascular-Safe modes. The HF generator automatically adjusts the output.

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Effect

The cutting or coagulation effect is preselected by the user in several intensities. The number of intensity levels depends on the mode selected (see table p. 59).

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Activation types

Three footswitch types are available to activate the HF energy. Single-pedal, dual-pedal and triple-pedal footswitches. Activation is assigned to an output socket. For unipolar sockets there is the option of using the fingerswitch of certain instruments for activation. This function is always activated and can only be switched off with the modes 'CUT OFF' and 'COAG OFF'. In certain bipolar modes the activation types 'AUTOSTART 1' and 'AUTOSTART 2' are also available. When contact is made with the tissue, depending on the impedance, the HF generator is automatically activated after a time delay defined by the user.

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Factory settings of the HF modes (procedure groups URO and GYN)

Predefined procedures for unit versions (High End) 125, U125

Name and description	Procedure group	Procedure no.	HF socket	Mode	Output [Watt]	Effect	Type of activation
Urology, Saline-Bi-TUR large loop	URO	100	Bi-MF	Saline-C-Cut	–	5	Double pedal
				Saline-Coag	–	6	
++Urology, Saline-Bi-TUR large loop		101*	Bi-MF	Saline-C-Cut++	–	5	Double pedal
				Saline-Coag++	–	6	
Urology C-Cut Unipolar-TUR-P		102	Unipolar	C-Cut	200	8	Double pedal
				Spray Coag	60	2	
Urology TOP-Cut Unipolar-TUR-P		103	Unipolar	TOP Cut	300	4	Double pedal
				Spray Coag	70	2	
Urology Unipolar-TUR-B		104	Unipolar	TOP Cut	120	4	Double pedal
				Spray Coag	60	2	
Urology – Circumcision	105	Bipolar	Bipolar Soft Coag	35	4	Double pedal	
Gynecology Unipolar-Hyst.	GYN	200	Unipolar	TOP Cut	140	4	Double pedal
				Forced Coag	50	3	
Gynecology – LASH - SupraLoop		201	Unipolar	TOP Cut	160	4	Double pedal
				Standard Coag	60	3	
		Bipolar	Bipolar Soft Coag	40	4	Single pedal	
Gynecology – Conization		202	Unipolar	TOP Cut	50	4	Double pedal
				Forced Coag	60	3	
Gynecology – Tube sterilization		203	Bipolar	Bipolar Soft Coag	40	4	Double pedal
Gynecology, Saline-Bi-Hyst. small loop		204	Bi-MF	Saline-C-Cut	–	5	Double pedal
				Saline-Coag	–	5	
Gynecology, Saline-Bi-Hyst. large loop		205	Bi-MF	Saline-C-Cut	–	5	Double pedal
				Saline-Coag	–	6	
Gynecology, Saline-Time-Hyst. small loop		206	Bi-MF	Saline-Time-C-Cut	–	5	Double pedal
				Saline-Time-Coag	–	6	
Gynecology Saline-Bi-Needle		207	Bipolar	Bipolar Cut	60	5	Double pedal
				Bipolar Soft Coag	40	4	
Gynecology Saline-TVE-Bi-Forceps		208	Bipolar	Bipolar Soft Coag	35	5	Double pedal

* This procedure is only available as a factory setting for the unit versions 205352 21x125 and 205352 23x125. The versions 205352 20x125 and 205352 22x125 can be updated for a surcharge.

Factory settings of the HF modes (procedure groups LAP, Gastro, open SUR, PED, ART, ENT, and VS)

Predefined procedures for unit versions (High End) 125, U125

Name and description	Procedure group	Procedure no.	HF socket	Mode	Output [Watt]	Effect	Type of activation
Laparoscopy RoBi & unipolar inst.	LAP	300	Bipolar	Bipolar Soft Coag	35	4	Single pedal
			Unipolar	TOP Cut	35	3	Double pedal
				Forced Coag	40	3	
Laparoscopy – Cholecystectomy		301	Bipolar	Bipolar Soft Coag	40	4	Single pedal
			Unipolar	Power Cut	60	3	Double pedal
				Forced Coag	50	3	
Laparoscopy – Hernia		302	Bipolar	Bipolar Soft Coag	40	4	Single pedal
			Unipolar	Power Cut	60	3	Double pedal
				Forced Coag	50	3	
Laparoscopy – Bi-Vascular-Safe + RoBi		303**	Bipolar	Bipolar Soft Coag	40	4	Single pedal
			Bi-MF	Bi-Vascular-Safe	–	2	Double pedal
Gastroscopy, unipolar polypectomy	Gastro	400	Unipolar	Gastro Cut	120	3	Double pedal
				Forced Coag	50	2	
Gastroscopy, unipolar papillotomy		401	Unipolar	Papillo Cut	100	2	Double pedal
				Standard Coag	40	2	
Open surgery Unipolar handle bipolar forceps	open SUR	500	Unipolar	TOP Cut	70	2	Double pedal
				Forced Coag	40	2	
			Bipolar	Bipolar Soft Coag	35	4	Single pedal
Open surgery, vessel sealing		501**	Bi-MF	Bi-Vascular-Safe	–	4	Single pedal
			Unipolar	TOP Cut	35	3	Double pedal
				Forced Coag	40	3	
Pediatrics – RoBi & unipolar Inst.	PED	600	Bipolar	Bipolar Soft Coag	30	3	Single pedal
			Unipolar	TOP Cut	30	3	Double pedal
				Forced Coag	35	3	
Arthroscopy – Hook electrode	ART	700	Unipolar	TOP Cut	50	3	Double pedal
				Standard Coag	45	3	
ENT, bipolar forceps	ENT	800	Bipolar	Bipolar Soft Coag	30	3	Double pedal
Visceral surgery	VS	900**	Unipolar	TOP Cut	50	3	–
				Forced Coag	40	2	
			Bi-MF	Bi-Vascular-Safe	–	4	Single pedal
			Bipolar	Bipolar Soft Coag	40	4	Double pedal

** This procedure is only available as a factory setting for the unit versions 205352 22x125 and 205352 23x125. The versions 205352 20x125 and 205352 21x125 can be updated for a surcharge.

Factory settings of the HF modes (procedure groups URO, GYN, LAP, Gastro, open SUR, PED, ART, and ENT)

Predefined procedures for unit versions (standard) 122, U122

Name and description	Procedure group	Procedure no.	HF socket	Mode	Output [Watt]	Effect	Type of activation
Urology C-Cut Unipolar-TUR-P	URO	100	Unipolar	C-Cut	200	8	Double pedal
				Spray Coag	60	2	
Urology TOP-Cut Unipolar-TUR-P		101	Unipolar	TOP Cut	300	4	Double pedal
				Spray Coag	70	2	
Urology Unipolar-TUR-B	URO	102	Unipolar	TOP Cut	120	4	Double pedal
				Spray Coag	60	2	
Urology – Circumcision		103	Bipolar 1	Bipolar Soft Coag	35	4	Double pedal
Gynecology – Unipolar HSK	GYN	200	Unipolar	TOP Cut	140	4	Double pedal
				Forced Coag	50	3	
Gynecology – LASH - SupraLoop		201	Unipolar	TOP Cut	160	4	Double pedal
				Standard Coag	60	3	
			Bipolar 1	Bipolar Soft Coag	40	4	Single pedal
Gynecology – Conisation		202	Unipolar	TOP Cut	50	4	Double pedal
				Forced Coag	60	3	
Gynecology – Tube sterilization		203	Bipolar 1	Bipolar Soft Coag	40	4	Double pedal
Gynecology Saline-Bi-Needle	GYN	204	Bipolar 1	Bipolar Cut	60	5	Double pedal
				Bipolar Soft Coag	40	4	
Gynecology Saline-TVE-Bi-Forceps		205	Bipolar 1	Bipolar Soft Coag	35	5	Double pedal
Laparoscopy RoBi & unipolar inst.	LAP	300	Bipolar 1	Bipolar Soft Coag	35	4	Single pedal
			Unipolar	TOP Cut	35	3	
				Forced Coag	40	3	Double pedal
Laparoscopy – Cholecystectomy		301	Bipolar 1	Bipolar Soft Coag	40	4	Single pedal
			Unipolar	Power Cut	60	3	
				Forced Coag	50	3	Double pedal
Laparoscopy – Hernia	LAP	302	Bipolar 1	Bipolar Soft Coag	40	4	Single pedal
			Unipolar	Power Cut	60	3	
				Forced Coag	50	3	Double pedal
Gastroscopy, unipolar polypectomy	Gastro	400	Unipolar	Gastro Cut	120	3	Double pedal
				Forced Coag	50	2	
Gastroscopy, unipolar papillotomy	Gastro	401	Unipolar	Papillo Cut	100	2	Double pedal
				Standard Coag	40	2	
				TOP Cut	70	2	
Open surgery - monopolar handle, bipolar forceps	open SUR	500	Unipolar	Forced Coag	40	2	Double pedal
			Bipolar 1	Bipolar Soft Coag	35	4	
			Bipolar 2	Bipolar Soft Coag	35	4	Single pedal
			Bipolar 1	Bipolar Soft Coag	30	3	Single pedal
Pediatrics – RoBi & unipolar Inst.	PED	600	Unipolar	TOP Cut	30	3	Double pedal
				Forced Coag	35	3	
Arthroscopy – hook electrode	ART	700	Unipolar	TOP Cut	50	3	Double pedal
				Standard Coag	45	3	
HNO, bipolar forceps	ENT	800	Bipolar	Bipolar Soft Coag	30	3	Double pedal

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Geeignete Elektroden

Nadelelektroden, Messerelektroden, Spatelektroden, Schlingenelektroden.

Cutting modes**POWER Cut***Properties*

Reproducible voltage-controlled atraumatic cutting, minimizing trauma to tissue, low to medium level of hemostasis.

The POWER Cut mode is equipped with PPS (Power Peak System). The AUTOCON® II 400 is equipped with an automatic power control system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arc necessary for the cutting quality selected even in the case of low-resistance loads. This is particularly important when the cutting wire electrode is pressed firmly against the tissue to be cut before activating the HF generator so that the cutting wire electrode covers a relatively large area and is therefore in low-resistance contact with the tissue. This is, for example, the rule with TUR and endoscopic polypectomy. In such cases the HF generator has to provide an above-average level of power on activation to ensure incision can be made without delay. Thanks to this device the average power can be limited to relatively low levels, thus reducing the risk of unintentional thermal damage to tissue.

Fields of application

All cutting procedures in tissue with good electric conductivity: e.g. muscle tissue, vascularized tissue. Dissection or cutting fine structures.

Suitable electrodes

Needle electrodes, knife electrodes, spatula electrodes, cutting loops.

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TOP Cut

Properties

Arc-controlled reproducible cutting which is atraumatic to tissue, in particular to tissue varying in nature and with poor electric conductivity.

PPS (Power Peak System)

The TOP Cut mode is equipped with PPS. The AUTOCON® II 400 is equipped with an automatic power control system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arc necessary for the cutting quality selected even in the case of low-resistance loads. This is particularly important when the cutting wire electrode is pressed firmly against the tissue to be cut before activating the HF generator so that the cutting wire electrode covers a relatively large area and is therefore in low-resistance contact with the tissue. This is, for example, the rule with TUR and endoscopic polypectomy. In such cases the HF generator has to provide an above-average level of power on activation to ensure incision can be made without delay. Thanks to this device the average power can be limited to relatively low levels, thus reducing the risk of unintentional thermal damage to tissue.

Fields of application

Includes the cutting of fatty structures and cutting under water, e.g. with TUR-P.

Suitable electrodes

Knife electrodes, spatula electrodes, cutting loops.

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C-Cut

Properties

High level of hemostasis with slightly arrested cutting.

Fields of application

E.g. cutting with low levels of bleeding during unipolar endoscopy with TUR-P and TUR-B, hysteroscopy.

Suitable electrodes

Electrodes with a large contact surface: knife and spatula electrodes, cutting loops and ribbon snare electrodes.

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PAPILLO Cut

Properties

Cutting consists of alternating cutting and coagulation phases. Cutting is easily controlled and is characterized by a reproducible and preselectable coagulation property that is specially suited to papillotomy during cutting.

Fields of application

Endoscopic interventions in which alternating cutting and coagulation is required with a single activation.

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SALINE-C-Cut

Properties

Saline-Time C-Cut is equipped with PPS system. Cutting current which is used with conductive NaCl irrigation solution in combination with bipolar KARL STORZ resectoscopes for cutting with low levels of bleeding during TUR-P and TUR-B and also surgical hysteroscopy with conductive NaCl irrigation solution.

PPS (Power Peak System)

Fields of application

E.g. cutting with low levels of bleeding during TUR-P and TUR-B and surgical hysteroscopy.

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SALINE-C-Cut ++

Features

Reproducible, atraumatic incisions. The effect intensities are used to set the degree of hemostasis at the incision edge.

PPS (Power Peak System)

Fields of application

Cutting procedures in bipolar resection.

PPS (Power Peak System)

The Saline C-Cut ++ mode is equipped with PPS. The initial cutting phase can pose a particular problem in cutting. This is especially true if the cutting electrodes have limited tissue contact when the HF generator is activated. This is the case in TUR, for example. In such cases, the HF generator has to provide an above-average level of power to ensure incision can be made without delay. The AUTOCON® II 400 is equipped with an automatic power regulation system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arcs necessary for the cutting quality selected even in the case of low-resistance loads. Thanks to this feature, the average output can be limited to relatively low levels, which corresponds to an improvement in safety against unintentional thermal tissue damage.

Suitable electrodes

Bipolar resectoscopes connected to the Bi-Saline resection socket on the resection auxiliary module with the bipolar cable 27176 LEB/LEBL.

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PPS (Power Peak System)

Fields of application

Cutting procedures in bipolar resection. PPS (Power Peak System) The Saline Time C-Cut ++ mode is equipped with PPS. During use, the Saline-Time-C-Cut ++ is emitted with automatic timed switch-off from 0.1 to 0.9 seconds. With premature disabling of the footpedal, the power output is immediately interrupted. During the activation time, the Saline-Time-C-Cut++ current is emitted when the yellow footpedal is depressed. Regulation of the power output when cutting with low tissue contact is the same as for Saline-C-Cut ++.

Suitable electrodes

Bipolar resectoscopes connected to the Bi-Saline Resection socket on the 'Resection' auxiliary module with the bipolar cable 27176 LEB/LEBL.

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BIPOLAR Cut

Properties

Cutting current which only flows directly around the distal end of the applicator. The effect intensities can be used to vary the degree of hemostasis at the incision edge.

The BIPOLAR CUT mode is equipped with PPS. The AUTOCON® II 400 is equipped with an automatic power control system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arc necessary for the cutting quality selected even in the case of low-resistance loads. This is particularly important when the cutting wire electrode is pressed firmly against the tissue to be cut before activating the HF generator so that the cutting wire electrode covers a relatively large area and is therefore in low-resistance contact with the tissue. This is, for example, the rule with TUR and endoscopic polypectomy. In such cases the HF generator has to provide an above-average level of power on activation to ensure incision can be made without delay. Thanks to this device the average power can be limited to relatively low levels, thus reducing the risk of unintentional thermal damage to tissue.

Suitable electrodes

Special applicators (bipolar electrodes with a rigid or height-adjustable cutting needle) during laparoscopy, neurosurgery and ENT.

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STANDARD Coag

Properties

High level of slough formation in tissue is prevented, adhesion of electrode to tissue is reduced. If you wish to make use of the potentially great 'coagulation depth' of STANDARD Coag, select a low effect intensity and coagulate for an extended period of time. If you can only coagulate for a brief period, select a high effect intensity.

Fields of application

For almost all operations which require fast and reliable coagulation or where adhesion of the electrode would have an adverse effect on the coagulation process.

Suitable electrodes

Electrodes with a large contact area, e.g. ball electrodes for deep coagulation.

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SPRAY Coag

Properties

No-contact efficient surface coagulation with increased voltage, low insertion depths. Automatic power output within preselected limits.

Fields of application

Coagulation of various types of bleeding.

Suitable electrodes

Knife electrodes, lancet-shaped electrodes, also cutting loops and ribbon snare electrodes with unipolar TUR.

WARNING: Only insulated unipolar metal forceps should be used with clamp coagulation.

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BIPOLAR SOFT Coag

Properties

Low voltages, slough formation of tissue is prevented, adhesion of electrode to tissue is greatly reduced. If you wish to make full use of the potentially great coagulation depth of BIPOLAR SOFT Coag, select a low effect intensity and coagulate for an extended period of time. If you can only coagulate for a brief period, select a high effect intensity. In comparison with other Coag modes you will still achieve a great coagulation depth, although you will not make full use of the potential coagulation depth of BIPOLAR SOFT Coag.

BIPOLAR SOFT Coag with AUTOSTART

In the window Activation type you can select the AUTOSTART function for BIPOLAR SOFT Coag. This additional function starts coagulation automatically after a defined time period when the instrument touches the tissue. Touching conductive metal material with the bare tips of the forceps, for example, will always initiate activation of the HF current with the additional AUTO START function!

WARNING: When introducing bipolar instruments into a metal trocar, the Auto Start function should be disabled beforehand!

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BIPOLAR SOFT Coag with AUTOSTOP

Properties

Identical to BIPOLAR SOFT Coag; however, it is firmly coupled with the AUTOSTOP function. AUTOSTOP ends activation automatically before the tissue adheres to the instrument or slough formation commences.

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SALINE-Coag

Properties

Coagulation mode for use with conductive NaCl irrigation solution in combination with KARL STORZ bipolar resectoscopes for hemostasis with TUR-P and TUR-B interventions as well as hysteroscopy. Here, the active resection loop must make contact with the bleeding vessel and the HF current activated at the same time.

Fields of application

E.g. for hemostasis with TUR-P and TUR-B interventions as well as hysteroscopy.

Suitable electrodes

Cutting loops recommended by KARL STORZ, vaporization electrodes, needle electrodes for bipolar resectoscopes from KARL STORZ with NaCl irrigation solution.

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Bi-VASCULAR-SAFE

i NOTE: Bi-Vascular-Safe is an optional mode that, if required, can be updated for a surcharge.

Features

Special COAG mode for the KARL STORZ Bi-Vascular-Clamp instruments (bipolar clamp for open surgery) and special RoBi® instruments with non-detachable MF cable. Using the four effect levels, you can dose the coagulation power exactly for the respective tissue type. The AUTO STOP function is calibrated to Bi-Vascular-Safe and automatically ends the activation when the best coagulation result is attained. The output is regulated completely automatically; no manual specifications for the power limitation are required.

Modulation

Bi-Vascular-Safe is a modulated form of current with alternating pulse and pause periods. The modulation levels can be changed under 'Mode settings'.

Fields of application

Open surgery, laparoscopy, and endoscopy.

Suitable electrodes

Please contact your KARL STORZ authorized dealers for further advice and information.

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Symbols for effects which can be set

Mode	Visualization in 'Parameter settings'								Visualization in 'Socket information'	Visualization in 'Mode selection'
	Effect 1	Effect 2	Effect 3	Effect 4	Effect 5	Effect 6	Effect 7	Effect 8		
Coag off										
TOP Cut										
POWER Cut & LAP C-Cut										
Bipolar Cut										
C-Cut										
Saline-C-Cut & Saline-Time-C-Cut & Saline-C-Cut++ & Saline-Time-C-Cut++										
Snare Cut										
Papillo Cut										

Mode	Visualization in 'Parameter settings'								Visualization in 'Socket information'	Visualization in 'Mode selection'
	Effect 1	Effect 2	Effect 3	Effect 4	Effect 5	Effect 6	Effect 7	Effect 8		
Standard Coag										
Forced Coag										
Spray Coag										
Bipolar Soft Coag										
Bi-Soft Autostop										
Bi-Vascular-Safe as of generator control software 2.00										
Bi-Vascular-Safe as of generator control software 20.00										
Saline-Coag & Saline-Time-Coag & Saline-Coag++ & Saline-Time-Coag++										

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























Unit versions:

For the AUTOCON® II 400 the following equipment versions are available with different socket assignments:

Application	Standard: Uni-/Bipolar “-” = 200-240V~	High End “-” = 200-240V~	Standard: Uni-/Bipolar “U” = 100-120V~	High End “U” = 100-120V~
Unit version	-122	-125	U122	U125
Article no. of device	20 5352 20-122	20 5352 20-125 Basic unit	20 5352 20U122	20 5352 20U125 Basic unit
		20 5352 21-125 Resection module		20 5352 21U125 Resection module
		20 5352 22-125 BiVascular		20 5352 22U125 BiVascular
		20 5352 23-125 Resection module and BiVascular		20 5352 23U125 Resection module and BiVascular
SET no. of device	20 5352 01-122	20 5352 01/02/03/04-125	20 5352 01U122	20 5352 01/02/03/04U125
Socket position				
1	Bipolar combo socket 	Bipolar combo socket 	Bipolar combo socket 	Bipolar combo socket 
2	Bipolar combo socket 	MF 	Bipolar combo socket 	MF 
3	Unipolar 3-pin and KS 	Unipolar 3-P and KS 	Unipolar 3-P and KS 	Unipolar 3-P and KS 
4	NE 2-pin and 6.3mm jack 	NE 2-pin and 6.3mm jack 	NE 2-pin and 6.3mm jack 	NE 2-pin and 6.3mm jack 

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Unpacking the equipment

Carefully remove the AUTOCON® II 400 and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Installing and connecting up the unit

! WARNING: The AUTOCON® II 400 HF surgery unit including the accessories connected may be used only in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations. This unit is not intended for use in hazardous zones. Do not operate the unit within the demarcated hazard zone while explosive anesthetic gases are in use.

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
Place the unit on a level surface. If installed on an unlevel surface, the unit must be secured against falling using suitable media. The unit might otherwise fall from the surface and injure people. This applies in particular to facilities that are mobile and/or can be swiveled such as ceiling supply units and equipment carts!

! WARNING: HF surgery equipment generates radiofrequency and electromagnetic fields that can interfere with the operation of sensitive electronic equipment. HF surgical equipment, and particularly its cabling, should thus never be installed close to any equipment sensitive to electromagnetic interference. Parallel cable routing with cables from other units should thus be avoided as a general rule, as HF cables emit HF emissions when activated, which may interfere with other equipment.

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
The instrument is equipped with a connector (8) for potential equilibration. The instrument's ground line should be installed by a qualified electrician.

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Checking selected line voltage (115/230 V)

NOTE: Before you open the line fuse holder, remove the power cord.

The AUTOCON® II 400 units are available in two voltages: 115 V and 230 V. You can tell the 115 V variants by the 'U' in the article number and they may only be operated at 115 V. The 230 V variants have a '-' in the article number on the manufacturer's identification plate and may only be operated at 230 V. How to check whether the fuse ratings and the voltage selection module are correct is described below.

If necessary, remove the cover of the line fuse holder/voltage selection module (see 'Fuse replacement' section). Take out the line fuse holder/voltage selection module. Turn through 180° until the indicator (115 or 230 V) matches the local voltage supply.

The setting 115 V is valid for a voltage range of 100–120 VAC, and the setting 230 V for a range of 220–240 VAC.

Make sure that the fuses used are suitable for the line voltage selected and are as stated on the manufacturer's identification plate:

205352 20-122/205352 2x-125; x = 0, 1, 2, 3
230 V: T4 AH250V (Art. no. 2027690)

205352 20U122/205352 2xU125; x = 0, 1, 2, 3
115 V: T8 AH250V (Art. no. 2028090)

Insert line fuse holder/voltage selection module back into unit, close cover.

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
Installation and operating instructions

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Connect power cord when the unit is switched off. Insert power cord into power cord receptacle ⑤ as far as it will go.

! WARNING: Make sure to use the correct line voltage and power cord type.

i NOTE: The instrument may only be operated with the power cord delivered by the manufacturer or a power cord of similar quality ("Hospital Grade Plug"), which has a national inspection seal.

Only operate the unit with the voltage stated on the identification plate.

! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

i NOTE: Please ensure that access to the power plug is ensured at the site of installation.

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

Insert the connecting cable of the footswitch into socket ⑤. Secure cable using screwed connection (union nut).

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SCB

i NOTE: To prevent the SCB connecting cable being pulled out accidentally, the SCB connector possesses a protection device. Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ①. Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control unit or other SCB units (see SCB control NEO System Instruction Manual).

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Operating the unit

Operating the unit for the first time

The AUTOCON® II 400 is delivered as follows:

Language: English

Neutral electrode: any

Maximum activation time: 99 sec.

Predefined procedures:

For every unit 100 user procedures can be saved. The memory locations 0 to 99 are reserved for this purpose. Every unit is equipped with several predefined procedures which cannot be permanently modified or deleted. They occupy memory locations as of 100.

ⓘ NOTE: The predefined procedures are divided into groups depending on medical specialties. The individual groups can be shown/hidden.

The predefined procedures set in the factory are divided into the following user groups:

- Urology
- Gynecology
- Laparoscopy
- Gastroenterology
- Open surgery
- Pediatrics
- Orthopedics and sports medicine
- ENT
- Visceral surgery

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Switch on the power switch.

When the unit is switched on, the logo screen is first displayed for approximately 5 seconds, then a white screen appears for approx. 2 seconds.

During power-up the unit performs a self-test. Here it checks important system components such as the HF generator, output sockets, touch screen, internal memory components, etc. and indicates any malfunction via an error message. The two signal paths for alarms, display and the audio system must be checked by the user whenever the unit is powered up.

NOTE: During power-up of the unit, the KARL STORZ logo must firstly be displayed followed by the white screen. In the event of any error the unit must not be used.

WARNING: When powering up the unit, pay attention to the acoustic availability signal.

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

The power-up message is followed by the 'Procedure List' screen. This lists the predefined user procedures and the procedures defined by the user himself. In routine operation the procedure information is called up here. If button ① is activated, the user will see the screen 'Device Setup: Select Option' (when the unit is operated for the first time, English is preset).

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Selection of language and user settings

The button 'User Settings' ②4 can be used to access the selection list of the user settings.

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

The arrow keys (2) next to the field 'Language' can be used to select the language. With this screen the user can also select the brightness of the display, the volume of the alarm system and the activation signals, the keyclick volume and the date format. If the user activates button (3), the properties selected will be permanently saved. If the user activates button (4), he will return to the 'User Settings' screen without saving the values. If the user continues to activate button (4), he will return to the 'Procedure List'.

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Selecting and modifying a predefined procedure

One of the predefined procedures can be selected on the 'Procedure List' screen by activating the relevant assignment buttons.

NOTE: The most recently used procedure is shown at the top of the selection list on the first page and is highlighted with a green dot in the top left corner.

The two levels are only differentiated by separate buttons with the supplement 'A' or 'B' after the procedure number if a level 'B' exists for a procedure number. A procedure level 'B' cannot be marked as 'last used', i.e., the order in the procedure list is always 'A' and then 'B', irrespective of which level was selected last (see also p. 31).

The procedures are ordered according to the procedure number in the order thereafter. The predefined procedures cannot be deleted, but are divided into medical specialties and can be selected for display. Their numbering begins in blocks of one-hundred (100, 200, 300...) respectively. In newly delivered units, it is only possible to select predefined procedures and all of the medical specialties are selected for display.

The button (31) can be used to open the 'Configuration of the procedure list' screen (see p. 81). This button is always shown after the last procedure button. If all the medical specialties are selected for display, this is the only button shown.

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Once a predefined procedure is selected on the 'Procedure list' screen, the 'Procedure information' screen is shown, in which all the parameters of all the sockets are visible. To amend the parameters the appropriate socket button must be activated e.g. the bipolar socket button ⑳.

Within the socket, the mode, the power and effect settings as well as the footswitch function can be set and changed.

- ① **NOTE:** With the procedure button ⑱ the user will return to the procedure list.
- ② **NOTE:** An MF socket in which no cable is inserted is identified as 'MF socket'. The identification changes to 'Bipolar MF socket' if a bipolar cable is connected and 'Monopolar MF socket' if a monopolar cable is connected.
- ③ **NOTE:** In the Saline modes and Bi-Vascular-Safe mode, no power limitation is specified, just the effect. Consequently, the corresponding value is not shown on the MF socket display.

Once a socket button is selected in the 'Procedure information' screen, the 'Socket information' screen appears with the corresponding parameters. To modify the operating parameter, activate the socket of the relevant operating mode, 'Cut' (yellow, ⑳) or 'Coag' (blue, ㉑).

- ④ **NOTE:** The buttons ⑤, ⑥ and ⑦ allow the user to switch between the individual sockets.
- ⑤ **NOTE:** If connections are made on more than one socket of the same type, then a running number (e.g., BI 1, BI 2) appears on the socket buttons alongside the type. An MF socket in which no cable is inserted is identified as 'MF'. The identification changes to 'BI MF' if a bipolar cable is connected and 'MONO MF' if a monopolar cable is connected. This also applies for the following menus.

If a socket is selected to which no activation type is currently assigned, the button with the symbol of a crossed-out footswitch appears.

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Once the desired operating mode button e.g., 'Cut' (yellow, ②) has been selected in the 'Socket information' screen, the 'Parameter settings' screen appears.

The parameters (effect or power limitation) for cutting or coagulating can be changed with the 'Plus/Minus' buttons. The newly selected values are active but not permanently saved. They can be permanently saved by following the instructions specified under the item 'Saving a personal self-defined procedure' (p. 73).

NOTE: The user can return to the 'Socket information' screen via button ④.

The operating mode button ② allows the user to return to the 'mode selection' screen.

NOTE: There is a maximum peak voltage value Upmax for every mode, depending on the effect. The respective current value is displayed directly above the pictogram (see arrow).

The limit values for effect and power limitation depend on the current mode (see Section 'Technical description' as of page 106).

In the Saline modes and Bi-Vascular-Safe mode, no power limitation is specified, just the effect. This means the corresponding range is not shown in the display as these modes are bipolar modes and not unipolar modes. Neutral electrode monitoring is not required in bipolar modes.

NOTE: During a Cut activation, the pictogram in the center of the image is shown in yellow; during a Coag activation, in blue. A separate pictogram is assigned to every effect intensity in this level. However, these differ only in the detail within a mode.

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
In the 'Mode selection' screen, the user can select the operating mode as required for the relevant socket. The different sockets offer their own selection for the relevant modes 'Cutting' and 'Coagulation'. Activating the button will switch the unit to Active mode. This is indicated by the green indicator in the left-hand corner of the button. The unit automatically switches back to the 'Socket Information' screen.

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
In the 'Socket information' screen, the buttons ⑤, ⑥ and ⑦ can be used to select the output parameters of the other sockets. The socket activated in each case is indicated by a green indicator in the top left corner of the button. At the same time, the pictogram showing the assigned activation type ⑩ is also lit up. The neutral electrode button ⑧ can be used to access the 'Neutral Electrode' screen.

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In the 'Neutral electrode' screen, the contact quality of an applied neutral electrode is specified in Ohm. An acceptable value is shown in green, an unacceptable value in red.

NOTE: If the neonatal monitoring is switched on, '(Neonatal)' is also added in a second line.

The button ④ is used to navigate back to the 'Socket Information' screen, and from there clicking ④ again takes you back to the 'Procedure Information' screen.

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
Saving a personal self-defined procedure


If parameters have been amended in a 'predefined procedure' as instructed under the steps 'Selecting and modifying a predefined procedure' (p. 69), they can be permanently saved for subsequent applications. The button 11 should be activated for this purpose.

How to overwrite and delete a procedure is described from p. 85 onwards.

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
The button 'Save as New Procedure' 15 can be used to access the menu to save the procedure.

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
The memory locations are displayed in ascending order. Memory locations that are already assigned are gray.

The Editor can be accessed by directly selecting the memory location or navigating with the arrow buttons and then selecting the memory location. The number of the memory location is shown on the right in the field for the button.

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
Firstly, a user name (first field) with at least one alphanumeric character must be entered, before you can move on to the procedure description (second field) with the TAB button. The TAB button remains gray until you enter a character. Once the user confirms the procedure by pressing button 17, it is saved to the specified memory location, and the screen switches to the procedure list.

NOTE: If upper-case letters have been activated, the green indicator on the shift button 17 lights up.

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
Saving a level 'B' with a self-defined procedure

With the optional triple-pedal footswitch 20017831 it is possible to switch backwards and forwards between two levels of a procedures of required. Here a self-defined procedure as described under 'Selecting and modifying a predefined procedure' (p. 69) should be created. This procedure must be changed as required by selecting new parameters and footswitch functions. When saving as described under 'Saving a personal self-defined procedure', the button 'Create Level B of Procedure' should be activated in the menu used to save a procedure. This allows the user to access the Editor for level 'B' of a procedure.

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
The two levels of the procedure will now appear one under the other in the procedure list.

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
Assigning activation type to a socket

In the 'Socket Information' screen of a procedure the relevant socket is assigned a footswitch activation. By adjusting or modifying a procedure in line with current requirements, the relevant footswitches, YELLOW for cutting and BLUE for coagulation, which are shown in the menu by colored pictograms next to the socket buttons, can be reassigned to another socket or blocked. To do so press button (10) next to the socket button. The button which can be activated has a light gray background and is assigned to the relevant socket. It is possible to connect two footswitches to the unit: a single-pedal footswitch and either a dual or triple-pedal footswitch. The single-pedal switch and thus also its pictogram only have one blue pedal. This means that it is only possible to activate the coagulation function.

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The footswitch activation selected, here yellow for cutting, is indicated by the green indicator (9). After selection the screen view automatically switches back to the 'Socket Information' screen, where footswitch activation can be assigned to additional sockets or blocked. With barred footswitch activation, HF energy can be released in unipolar mode via the finger-operated buttons on the instrument. For bipolar applications automatic activation can be selected. This is time-dependent and can be specified in the 'Device Setup' screen.

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Autostart function

NOTE: In the Autostart mode the output power limit is max. 50 W.

The Autostart function permits automatic activation of the coagulation current in the coagulation modes of the bipolar socket as soon as both branches of the coagulation electrode touch the tissue and if there is a defined impedance range. The Autostart function 1 allows coagulation to be activated from between 0.1 and 1.9 seconds in steps of 0.1 second. The Autostart function 2 provides for a time delay of between 2 and 10.0 seconds in steps of 0.5 seconds.


Autostart 1 is factory-set to a value of 1.0 seconds, and Autostart 2 to a value of 2.0 seconds. The user is free to reselect the time-to-activate delays at his own discretion.

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- Bipolar Soft Coag
- Bipolar Soft Coag mit Auto-Stop.

NOTE: The Autostart function can be used for the operation modes

- Bipolar Soft Coag
- Bipolar Soft Coag with Auto-Stop.

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The field for selection of the activation types is now shown next to the pictographic HF socket.

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The autostart buttons 1 and 2 with the pre-set time delays are now available for selection in the selection list next to the footpedal activator.

The permitted time selected under 'Application Settings' (see p. 79) is then shown.

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
! WARNING: The HF generator is now automatically activated when the tips of the forceps come into contact with tissue!

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
Time setting for bipolar Autostart activation

To set the time for bipolar Autostart activation, the button ① 'Setup Device' should be activated from the 'Procedure List' screen.

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
The button 'Application Settings' ⑯ can be used to switch to the 'Application Settings' screen for changing settings as necessary.

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
The required times can be set with the arrow buttons next to the fields Autostart 1 and Autostart 2. After selection, button ③ is used to save the parameters, and the system automatically returns to the 'Device Setup screen: Select Option'. With button ④ the user can move back without saving changes.

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
Selecting the neutral electrode

The type of neutral electrode is selected for all unipolar procedures of the AUTOCON® II 400. It can be adjusted as follows.

From the start screen 'Procedure List' activate the button ① 'Setup Device'.

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
The button 'Application Settings' ⑮ can be used to switch to the 'Application Settings' screen for changing settings as necessary.

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The required type can be selected with the arrow buttons next to the field 'NE Type'. The settings single-faced, double-faced/dynamic and 'any' can be selected for the connection of single-faced and double-faced neutral electrodes (for explanation, see pages 24-25, 27).

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The measuring range is from 0 - 199 ohms (lower and upper ends of bar). Greater resistances are shown in the numerical display by '----'.

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
① **NOTE:** If the neonatal monitoring is switched on, '(Neonatal)' is also added in a second line.

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Setting the maximum activation time

The maximum activation time can be adjusted according to the user's requirements. From the start screen 'Procedure List' activate the button ① 'Setup Device'. This setting is protected with a user password.


WARNING: Any change to the maximum activation time must be notified to all users.

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
The user can use the arrow buttons ② next to the field 'Max. Activation Time' to set the required time. When the arrow buttons ② are activated, the first screen to appear is that for entering the user password.

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NOTE: On delivery the user password is '12345'. It is not identical to the service password and can be changed as described in the section 'Changing the user password' (p. 83).

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
Entry of the password is secret, i.e. only asterisks appear in the display in place of the characters entered. After entering the valid password the user can confirm using button (3) and return to the 'Application Settings' screen where the buttons for selection of the maximum activation time have been enabled. After specifying the activation time this value can be permanently saved with button (3). To make any change to the activation time the valid password must be reentered.

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Connecting instruments

Depending on the type of unit, the unipolar and/or bipolar instruments can be connected to the sockets on the front using the appropriate cable.

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Connecting the footswitch

At the rear of the unit the user can connect the following as required:

- a single-pedal footswitch,
- a dual-pedal footswitch,
- a triple-pedal footswitch or
- a combination of a dual-pedal footswitch and a single-pedal footswitch and
- a combination of a triple-pedal footswitch and a single-pedal footswitch.


The 'Assigning activation type to a socket' is described on p. 75.

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
Activating HF energy

It is only possible to activate HF energy when the power and effect parameters are visible to the user. This is only possible from the 'Procedure Information', 'Socket Information' and 'Parameter Settings' screens.

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
Alarm and information messages

The AUTOCON® II 400 is equipped with an error control system. User errors and unit errors are displayed flashing (13). The codes are listed in the table on page 97.

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Further setting options:

Configuration of the procedure list

The predefined procedures are divided into groups of medical specialties. The user may select those he wishes to have displayed in the 'Procedure list' screen. The buttons have a toggle function. The green 'LED' signalizes that the corresponding specialty is selected for display.

NOTE: The 'Visceral Surgery' selection button is only available on the unit versions 20535222x125 and 20535223x125, i.e., following a successful upgrade to the Bi-Vascular-Safe mode (see Unit versions p. 62).

If all the specialties are deselected (and thus not displayed), only the button (31) 'Configuration of the procedure list' is shown in the procedure list. This button returns you to the screen.

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
Setting date and time

The AUTOCON® II 400 is supplied with the date and time set to 'Central European Time' CET.

From the start screen 'Procedure List' activate the button (1) Setup Device.

NOTE: The ticking clock shows that the screen is being updated. If the clock has stopped, the screen is frozen, and this error must be rectified before using the HF unit again.

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

With the button 23 'Date and Time' the user can switch to the screen used to set the date and time.

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
The settings can be made using the arrow buttons under the fields. After selection, button 3 is used to save the specifications, and the system automatically returns to the 'Device Setup: Select Option' screen.

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NOTE: The day of the week is calculated using the calendar. Selection is not necessary here.

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
The format for display of the date can also be selected. With the button 'User Settings' 24 the user can access the 'User settings' screen.

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
If the user activates button 3, the properties selected will be permanently stored. If the user activates button 4, he will return to the screen 'User Settings'. If the user activates this button again, he will return to the 'Procedure List'.

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
Changing the user password

If button ① is pressed from the menu 'Procedure List' the user can access the screen 'Device Setup: Select Option'.

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
The button 'Change User Password' ②5 can be used to access the 'Enter/Edit Password' screen.

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
Entry of the password is secret for security reasons – 'asterisks' appear in the display in place of the characters entered. The password has been factory-set to '12345'. The entry of a valid password takes the user to another entry menu, which only differs from the previous menu by having two entry lines. The text of the header changes to 'Input New Password'.

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
NOTE: The user password is not identical to the service password.

After the new password has been entered in line 1, pressing the button ⑫ will move the cursor to the second input line.

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
The text in the header changes to 'Confirm New Password'. After confirming the password, pressing button ③ will take the user back to 'Device Setup: Select Option' if the entry is valid, and the new password will be accepted. From here you can return to the 'Procedure List' screen with button ④. Activation of the button ③ with an invalid password will result in the message 'Invalid password'. The content of the two text lines is then deleted and the cursor is again positioned in the left-hand input line. Pressing the Back button ④ will take the user back to the Set-Up selection list, with the password remaining unchanged.

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Overwriting a procedure


Here there is the option of overwriting a procedure if the parameters of a procedure have not proved to be effective and are to be permanently amended.

The procedure to be overwritten is selected from the procedure list by pressing the appropriate button. This accesses the 'Procedures Information' screen.

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
Parameters are selected as described under 'Selecting and modifying a predefined procedure' (p. 69). From the overview of procedures, button ⑪ can be used to access the 'Storage Options' screen.

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
The button 'Overwrite Procedure' is used to determine the procedure displayed in the header to be overwritten. The procedure is overwritten using button ③, and the user is then taken to the 'Procedure List' screen.

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
Deleting a procedure

The procedure to be deleted is selected from the procedure list by pressing the appropriate button. This accesses the 'Procedure Information' screen.

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
From the 'Procedure Information' screen button 11 can be used to access the 'Storage options' screen.

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
The button 'Delete Procedure' 26 is used to determine the procedure displayed in the header to be deleted. The procedure is deleted using button 3, and the user is then taken to the procedure list.

User and procedure names can be changed using the option 'Rename Procedure' 27.

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Display of system events

The system events are listed in the logbook. From the 'Procedure List' screen, button 1 can be used to access the 'Device Setup: Select Option' screen.

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

The button 'System Log' ⑳ is used to access the logbook in which the system events are stored.

① NOTE: The last 200 system events are stored. A new entry overwrites the oldest entry. The entries in the logbook are retained permanent-ly even after the device has been switched off.

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
The system events can be viewed using the arrow buttons on the right of the Fig. With button ④ the user can return to the 'Device Setup Select: Option' screen.

From here you can return to the 'Procedure List' screen using button ④.

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
Display of device information

The device information states the serial number and software versions. From the menu 'Procedure List', button ① can be activated to access the menu 'Device Setup Select Option'.

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
The button 'Device Configuration' ㉔ is used to access the 'Device Information' screen, where the configuration data is stored.

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This screen displays the software status, the serial number, the expired operating hours, and the date of the next safety check, among other things. Optional modes, which are enabled, are also displayed here.


With button ④ the user can return to the 'Device Setup: Select Option' screen.

From here you can return to the 'Procedure List' screen using button ④.

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Mode specific settings


The 'Mode settings' can be accessed in the 'Device Setup: Select Option' screen via the 'Mode settings' button. In this screen it is possible to select the activation times (0.1 s - 1 s, in 100 ms steps) for the modes Saline-Time-Cut/Saline-Time-Cut++ and Saline-Time-Coag/Saline-Time-Coag++ using the UP/DOWN buttons. Activation is stopped according to the time selected regardless of whether the footswitch or fingerswitch is activated again.

In addition, the number of cutting impulses (1 – 5) can also be set separately for the gastro and papillo cut.

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The 'Down' button can be used to access further setting options:


The cutting speed (1 – 20) can be set separately for the gastro and papillo cut.

The modulation (5 – 300) can be set for the Bi-Vascular-Safe mode.

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
Service Programs

The button 'Service Programs' ③⑩ can be used to access the password-protected service area.

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Fuse replacement


Switch off the unit and remove the power cord from the electrical outlet. Only use fuses suitable for the voltage range stated on the manufacturer's identification plate.

! WARNING: Use only fuses of the correct rating. Replace defective fuses (see table) and reinsert line fuse holder/voltage selection module.

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Release cover from the line fuse holder/voltage selection module with a screwdriver or other suitable tool. Remove the line fuse holder/voltage selection module.

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! WARNING: Use only fuses of the correct rating.


Replace defective fuses (see table) and reinsert line fuse holder/voltage selection module.

! WARNING: Ensure that the display in the voltage selection module corresponds to the local power supply (115 oder 230 V).

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Reconnect the power cord. Test the unit for proper operation.

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Reprocessing

! WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

! WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

! WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

! WARNING: National laws and regulations must be observed.

i NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

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AUTOCON® II 400 and footswitch

! WARNING: Before any cleaning work is performed, the device is to be disconnected from the power source!

! CAUTION: Make absolutely sure that no liquid can get into the device housing.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

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Reprocessing of accessories

Extensive accompanying documentation is available for the accessories. It is essential that you refer to these documents for further information on handling and reprocessing.

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- Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.
- Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.
- Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, No. 1, pp. 23-29, 1995.

References for Cleaning and Sterilization

- Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.
- Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.
- Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, No. 1, pp. 23-29, 1995.

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Maintenance and safety check

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Safety check

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, safety checks/repeat inspections according to IEC 62353 must be carried out and recorded at least once a year by qualified personnel.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Protective ground resistance measured according to IEC 62353: For limit values, please refer to the current standard.
- Earth leakage current measured according to IEC 62353: For limit values, please refer to the current standard.
- Touch current measured according to IEC 62353: For limit values, please refer to the current standard.
- Patient leakage current measured according to IEC 62353: For limit values, please refer to the current standard.

Test for proper operation

The test for proper operation must include the following tests at least:

- Testing for proper operation of all visual and acoustic signals

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Maintenance

- Testing of the monitoring devices
- Testing of the function modes
- Testing of the automatic start mode
- Checking the activation of the HF power using the dual-pedal footswitch
- Testing the HF power in the various operating modes

Documentation

The safety check and the results must be recorded in the medical device logbook.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during this safety check/repeat inspection, the unit must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety checks, please refer to the current version of the relevant Service Manual.

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Servicing and repair

WARNING: Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

To bridge the repair period, you will, as a rule, receive a unit on loan which must then be returned to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

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In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

Manufacturer's warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**AUTOCON® II 400 Alarm and
information messages**

Code	Description	Rectification / explanation	Priority	Alarm conditions
A-0B	NE contact: Resistance too low	Check neutral electrode and reposition if necessary	Medium	physiological
A-0C	NE contact: Resistance too high	Check neutral electrode and reposition if necessary	Medium	physiological
A-10	End activation (AutoStop reached)	End activation of HF energy	Low	other
A-11	No activation signal	Have unit serviced.	High	technical
A-12	Activation during power up not permitted	Check footswitch, handle and cable	Medium	other
A-17	Simultaneous activation not permitted	Activate either via footswitch or handle	Medium	other
A-1C	Max. activation time exceeded	Cancel activation	Low	other
A-1D	Touch screen short circuit	Touch screen short circuit or pressed too long	Low	technical
A-1E	Touch screen pressed during power up	Touch screen must not be pressed during power-up	Low	technical
A-1F	NE symmetry critical (cancel activation)	Check neutral electrode and reposition if necessary	High	technical
A-22	AutoStop reached	Break contact with tissue	Low	other
A-85	Dual-pedal footswitch disconnected	Reconnect dual-pedal footswitch and check cable	Low	technical
A-89	Single-pedal footswitch disconnected	Reconnect single-pedal footswitch and check cable	Low	technical
A-94	Triple-pedal footswitch disconnected	Reconnect triple-pedal footswitch and check cable	Low	technical
A-95	Instrument detected at MF socket	New MF instrument at socket	Low	technical
A-96	Number of applications approx. nn	Follow-up message to A-95	Low	technical
A-9A	Incompatible MF instrument recognized or Bi-Vascular-Safe upgrade not installed	Only use compatible MF instruments Install Bi-Vascular-Safe upgrade	Low	technical
A-9F	MF instrument disconnected	MF instrument removed from socket	Low	technical
A-A3	Footswitch not assigned to a socket	Assign footswitch to a socket (p. 75)	Low	other
A-A4	2 footswitches of same type connected	Either connect one footswitch or two of a different type	Low	other
A-A9	First select procedure before activation	Select the required procedure (p. 69)	Low	other
A-AA	No mode assigned	Assign a mode to the procedure (p. 72)	Low	other
A-AB	No/Unsuitable MF instrument detected	Attempt to activate the socket where there was no or a defective MF instrument detected.	Medium	technical
A-AC	Autostart on tissue contact not possible	Release tissue contact, separate instrument tips, if necessary, assign Autostart	Medium	other
A-B0	NE symmetry not optimal	Check neutral electrode and reposition if necessary	Low	technical
A-B2	NE current limit exceeded (>300 mA)	The current across the neutral electrode exceeds 300 mA (risk of burns, particularly in neonatal applications)		
A-B7	Autostart only permitted up to 50W output	The max. output with Autostart should only be 50 W	Low	technical
A-BB	SC due	Arrange date for safety check	Medium	technical
C-84	Short-circuit between AE and NE	No impedance modification of tissue during coagulation	Low	technical
all others		Inform Service		

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POWER Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	500 ohms
Max. HF peak voltage	740 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	10 W to 300 W in 1 W steps
Max. power output at rated load resistance	300 W $\pm 20\%$

TOP Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	500 ohms
Max. HF peak voltage	1040 Vp (with electric arc)
Number of effects	8
Constancy of effects	automatic control of arc intensity
HF power limitation	10 W to 300 W in 1 W steps
Max. power output at rated load resistance	300 W $\pm 20\%$

LAP-C-Cut	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	
Effect 1-4:	3.2 (at $R_L = 500 \text{ ohms}$)
Effect 5+6:	3.3 (at $R_L = 500 \text{ ohms}$)
Effect 7+8:	3.6 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	500 ohms
Max. HF peak voltage	1450 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	10 W to 200 W in 1 W steps
Max. power output at rated load resistance	200 W $\pm 20\%$

PAPILLO Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz $\pm$ 10 %
Crest factor	1.4
Rated load resistance	75 ohms
Max. HF peak voltage	800 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
High frequency power limitation	10 W to 200 W in 1 W steps
Max. power output	200 W $\pm$ 20 %

GASTRO Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz $\pm$ 10 %
Crest factor	75 ohms
Rated load resistance	1.4
Max. HF peak voltage	880 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
High frequency power limitation	10 W to 200 W in 1 W steps
Max. power output	200 W $\pm$ 20 %

BIPOLAR Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500$ ohms) $\pm$ 10 %
Crest factor	1.4 (at $R_L = 500$ ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	740 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
High frequency power limitation	10 W to 100 W in 1 W steps
Max. power output	100 W $\pm$ 20 %

C-Cut	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	
Effect 1-4:	3.2 (at $R_L = 500 \text{ ohms}$)
Effect 5+6:	3.3 (at $R_L = 500 \text{ ohms}$)
Effect 7+8:	3.6 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	500 ohms
Max. HF peak voltage	1450 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	10 W to 200 W in 1 W steps
Max. power output at rated load resistance	200 W $\pm 20\%$

SALINE-C-Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	500 ohms
Max. HF peak voltage	770 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
Max. power output	370 W $+8\%/-20\%$

SALINE-C-Cut++	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	75 ohms
Max. HF peak voltage	490 Vp
Number of effects	8
Constancy of effects	automatic control of HF generator voltage supply
Max. power output at rated load resistance	300 Watt $\pm 20\%$

SALINE-TIME-C-Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	500 ohms
Max. HF peak voltage	770 Vp
Number of effects	8 (time 0.1 - 1 sec.)
Constancy of effects	automatic control of HF peak voltage
Max. power output	370 W $+8\%/-20\%$

SALINE-TIME-C-Cut++	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	75 ohms
Max. HF peak voltage	490 Vp
Number of effects	8 (time 0.1 - 1 sec.)
Constancy of effects	automatic control of HF generator voltage supply
Max. power output at rated load resistance	300 W $\pm 20\%$

STANDARD Coag	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	50 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	5 W to 200 W in 1 W steps
Max. power output at rated load resistance	200 W $\pm 20\%$

FORCED Coag	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500$ ohms) $\pm 10\%$
Crest factor	5.0 (at $R_L = 500$ ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	1800 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
HF power limitation	5 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W $\pm 20\%$

SPRAY Coag	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500$ ohms) $\pm 10\%$
Crest factor	7.4 (at $R_L = 500$ ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	4300 Vp
Number of effects	2
Constancy of effects	limitation of HF peak voltage
HF power limitation	5 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W $\pm 20\%$

BIPOLAR SOFT Coag / BIPOLAR SOFT with AUTOSTOP	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency 350 kHz	350 kHz (at $R_L = 500$ ohms) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500$ ohms)
Rated load resistance	75 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	1 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W $\pm 20\%$

BI-VASCULAR-Safe	
HF voltage form	modulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	25 ohms
Max. HF peak voltage	220 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
Max. HF power	300 W $\pm 20\%$

SALINE-Coag	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	75 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
Max. power output at rated load resistance	200 W $\pm 20\%$

SALINE-Coag++	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500 \text{ ohms}$) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500 \text{ ohms}$)
Rated load resistance	50 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF generator voltage supply
Max. power output at rated load resistance	200 W $\pm 20\%$

SALINE-TIME-Coag	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500$ ohms) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500$ ohms)
Rated load resistance	75 ohms
Max. HF peak voltage	190 Vp
Number of effects	8 (time 0.1 - 1 sec.)
Constancy of effects	automatic control of HF peak voltage
Max. power output at rated load resistance	200 W $\pm 20\%$

SALINE-TIME-Coag++	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at $R_L = 500$ ohms) $\pm 10\%$
Crest factor	1.4 (at $R_L = 500$ ohms)
Rated load resistance	50 ohms
Max. HF peak voltage	190 Vp
Number of effects	8 (time 0.1 - 1 sec.)
Constancy of effects	automatic control of HF generator voltage supply
Max. power output at rated load resistance	200 W $\pm 20\%$

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF

Directive compliance

(for 20535220-122/2053522x-125)

This medical device belongs to Class IIb

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

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NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Recommended accessories

Item	Cat. no.
BiVascularClamp 200 , total length 20 cm, with thermoisolation, for the open surgery, integrated MF cable length 4 m, Vp max 250V, for use with AUTOCON® II 400 in the versions 20 5352 22U125 / 20 5352 23U125 module and 20 5352 22-125 / 20 5352 23-125	20 5400 21
BiVascularClamp 210 , total length 21 cm, with wavy jaws, for the open surgery, integrated MF cable length 4 m, Vp max 250V, for use with AUTOCON® II 400 in the versions 20 5352 22U125 / 20 5352 23U125 module and 20 5352 22-125 / 20 5352 23-125	20 5400 22
BiVascularClamp 270 , total length 27 cm, with thermoisolation, for the open surgery, integrated MF cable length 4 m, Vp max 250V, for use with AUTOCON® II 400 in the versions 20 5352 22U125 / 20 5352 23U125 module and 20 5352 22-125 / 20 5352 23-125	20 5400 23

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Recommended accessories

Item	Cat. no.
Triple-pedal footswitch , for activation of coagulation and cutting function, with function for switching between mode levels A and B	20017831
Dual-pedal footswitch , for activation of coagulation and cutting function	20017830
Single-pedal footswitch , for activation of coagulation function	20017832
Neutral electrode , made of conductive silicone with 2 rubber straps for fastening, contact area A = 500 cm ² , for use with connecting cable 27806	27805
Connecting cable , to connect neutral electrodes 27805 and 860021 E, length 300 cm	27806
Neutral electrode , disposable, contact area subdivided, A = 169 cm ² , (pack of 50), connecting cable 27801 also required	27802
Connecting cable , to connect neutral electrode, disposable 27802 and 860021 S, length 500 cm	27801
Electrode handle , with 2 buttons to activate unipolar generator, yellow button: unipolar cutting, blue button: unipolar coagulation (cable 26 5200 45 also required)	26 5200 43
High frequency cable , for electrode handle 26 5200 43, length 500 cm	26 5200 45
Electrode handle , without buttons, with built-in connecting cable, length 300 cm. The footswitch is used to activate the high frequency surgical unit	26 5200 46
External hardware upgrade AUTOCON® II 400 , suitable for the unit versions -125	27176B
Upgrade for Bi-Vascular-Safe	- - -

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and operation.

The AUTOCON® II 400 model 205352 20* corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

* applies to the AUTOCON® II 400 model versions 205352 20-122, 205352 20U122, 205352 20-125, 205352 20U125, 205352 21-125, 205352 21U125, 205352 22-125, 205352 22U125, 205352 23-125, and 205352 23U125.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The AUTOCON® II 400 model 205352 20 is a Group 1 unit (as per CISPR 11).

NOTE: The tables and guidelines that are included in this Appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The AUTOCON® II 400 model 20535220 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the AUTOCON® II 400 model 20535220. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

NOTE: When activated, the AUTOCON® II 400 HF surgical unit may interfere with other devices in the immediate vicinity. Telling signs are e.g. image artifacts in the case of imaging devices or unusual fluctuations in measurement displays.

Such interference by an activated AUTOCON® II 400 HF surgical unit can be reduced by taking the measures mentioned in the note on p. 132.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
Handle BOWA 214-03	No		No	Handle
20530076	No		No	Neutral electrode
Handle MTP	No		No	Handle
Erbe 20489-100	Yes	4.5 m	No	Double pedal footswitch
Erbe 20188-100	Yes	4.5 m	No	One pedal footswitch
20530045	No	5.0 m	No	HF Cord
20530070	No	3.0 m	No	NE Cord
20090170	Yes	1.0 m	No	SCB Cord
PE	No	3.0 m	No	Potential equalization
Power Cord	No	3.0 m	No	Power Supply

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
AUTOCON® II 400 model 20535220 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® II 400 model 20535220 should make sure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	While in stand-by the AUTOCON® II 400 model 20535220 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The AUTOCON® II 400 model 20535220 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
AUTOCON® II 400 model 20535220 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® II 400 model 20535220 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T^* (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T^* (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
AUTOCON® II 400 model 20535220 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® II 400 model 20535220 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V_{rms} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V_{rms} 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the AUTOCON® II 400 model 20535220, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5 / 3] \sqrt{P}$ $d = [3.5 \cdot 3] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = [7 / 3] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 6

**Recommended separation distances between portable and mobile RF communications equipment and the
AUTOCON® II 400 model 20 5352 20**

The AUTOCON® II 400 model 20 5352 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseas.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufnategui 627 6º Piso
B1638 CAA Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytøften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info@karlstorz.es

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopiija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoszkóp Magyarorszag Kft.
Toberek utca 2. fsz. 17/b
HU-1112 Budapest, HUNGARY
Tel.: +36 18775070, Fax: +36 18775070
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece E.P.E.*
Patriarhou Grigorion E' 34
54639 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
*Repair & Service Subsidiary

*Repair & Service Subsidiary
KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC "Arman", office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Onon Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.

Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.

Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, No. 1, pp. 23-29, 1995.

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Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.

Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities. A Guide for Device Manufacturers. AAMI TIR № 12-1994.

Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, № 1, pp. 23-29, 1995.

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KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kazakhstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@itte.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 - 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty. Ltd.
174 Parramatta Rd.
Camperdown, N.S.W. 2050, Australia
Phone: +61 02 8594 9150, Fax: +61 02 8594 9199
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F, Chuanxing Mansory,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
10-18 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6FL, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

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