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
INTENDED USE

Urea nitrogen (BUN) liquid reagent is used for the kinetic quantitative determination of Urea Nitrogen (BUN) in human serum used in routine examination and monitoring of therapy and relapses.

SUMMARY AND PRINCIPLE

Urea is the principle waste product of protein catabolism. It is synthesized in the liver from ammonia which is produced as a result of the deamination of amino acids. Normally, urea nitrogen in the blood comprises only about 45% of the non-protein nitrogen. The importance of urea nitrogen determination is its value as an indicator of liver and kidney functions. Decreases in Blood Urea Nitrogen (BUN) are seen with nephritis, acute liver destruction, amyloidosis and pregnancy. Increases in BUN are encountered with acute and chronic nephritis, intestinal and urinary obstruction, uremia, metallic poisoning, pneumonia, Addison's disease, peritonitis, surgical shock and cardiac failure. This procedure is a modification of the method described by Sampson. Urea is catalytically converted to ammonium carbonate by the use of urease. The reaction rate is dependent upon the concentration of the influence of glutamic dehydrogenase. The rate of this second reaction is dependent upon the first and can be measured by the rate of conversion of NADH to NAD by the change of absorbency at 340 nm.

Urease

REAGENTS

BUN Liquid Reagents 1 and 2 and BUN Standard come in separate containers, and all reagents are clear, colorless liquid in ready to use format.

1. After combining BUN Liquid R1 and BUN Liquid R2 as directed the reagent contains:

TRIS Buffer, pH 7.8	100 mmol/L
2-Oxoglutarate	5 mmol/L
ADP	0.6 mmol/L
Urease	>20,000 U/L
GLDH	>1,500 U/L
NADH	0.25 mmol/L
Stabilizers, Preservatives	

2. Urea Nitrogen Standard (20 mg/dL)

WARNINGS AND PRECAUTIONS

- The reagents are for *in vitro* diagnostic use and are intended for professional use only.
- Normal precautions exercised in handling laboratory reagents should be followed.
- The reagents contain sodium azide, which may be toxic if ingested. Sodium azide may also react with lead and copper plumbing to form highly explosive metal azides.
- Refer to Material Safety Data Sheet for any updated risk, hazard, or safety information.

REAGENT PREPARATION

- The working reagent is prepared by mixing five (5) volumes of R1 with one (1) volume of R2 in a disposable container.

Example: 25 ml R1 + 5 ml R2

- Urea Nitrogen Standard is ready to use

REAGENT STORAGE AND STABILITY

Both the BUN reagents and standard are stable until the expiration date on their respective labels, when properly stored at 2-8°C and protected from light. Reagents should appear clear and colorless. Discard if either appears cloudy or contains particulate matter. The working reagent is stable for 2 weeks at 2-8°C. Protect the reagent from direct light. The working reagent should be discarded if the initial absorbance, read against distilled water at 340 nm, is below 1.000.

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- Spectrophotometer capable of absorbance reading at 340 nm
- Constant temperature block or bath, 37°C, or temperature controlled cuvette
- Accurate pipetting devices

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4. Test tubes
5. Interval timer

SPECIMEN COLLECTION AND STORAGE Non-hemolyzed serum is the specimen of choice. Whenever possible specimens should be separated and analyzed on the day of collection. Anticoagulants containing ammonium or fluoride salts must be avoided. Urea in serum is stable for up to 24 hours at room temperature (15-25°C), several days refrigerated at 2-8°C and for at least 2-3 months when frozen (-20°C).

INTERFERING SUBSTANCES

Fluoride and ammonia cause interference with the BUN assay. Blood collected in tubes containing ammonium heparinate should NOT be used. Young, et al⁴ provide a list of drugs and other substances that interfere with the determination of BUN.

MANUAL PROCEDURE

1. Prepare BUN working reagent according to instructions.
2. Pipette 1.0 mL of working reagent into tubes labeled "standard", "control", "patient", etc.
3. Pre-incubate all tubes at 37°C for at least five minutes.
4. Zero spectrophotometer at 340 nm with distilled water.
5. Add 10 µL (0.010 mL) of sample, mix and return to a thermo cuvette.
6. After exactly 30 seconds, read and record absorbance (A₁).
7. At exactly 60 seconds after reading (A₁), read and record absorbance (A₂).
8. Repeat steps 5-7 for all test specimens and standard.
9. Calculate change in absorbance (ΔA = A₂ - A₁) per minute.
10. See "Calculations".

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* HTI MULTI-PURPOSE CALIBRATOR MAY BE USED IN PLACE OF STANDARD.

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QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established BUN values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedural errors.

CALIBRATION

Calibration is required. The instrument manufacturer's calibration guidelines should be followed to calibrate your analyzer.

CALCULATIONS

Values are derived by comparing the absorbance change (ΔA) of the unknown (u) with that of a standard (s) identically treated.

$$\text{Serum BUN (mg/dL)} = \frac{\Delta A_u \text{ Unknown}}{\Delta A_s \text{ Std. or Cal.}} \times \text{Conc. Std.}$$

Where: A_u and A_s are the absorbance changes (decrease) of unknown and standard, respectively, and Conc. Std. is the concentration of standard (mg/dL).

Example: ΔA_u = 0.035, ΔA_s = 0.045, Conc. Std. = 20 mg/dL

$$\text{Serum BUN (mg/dL)} = \frac{0.035}{0.045} \times 20 = 15.6 \text{ mg/dL}$$

0.045

NOTE: To convert the results into SI units (mmol/L), multiply the result (mg/dL) by 0.357

LIMITATIONS

If the BUN value exceeds 80 mg/dL, the specimen should be diluted 2-fold (1+1) with distilled water; the assay repeated and results multiplied by the dilution factor of 2. BUN values for neonatal patients have not been established with this procedure.

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Urea (mg/dL) = BUN (mg/dL) × 2.14

Urea (mmol/L) = Urea (mg/dL) × 0.167

This range should serve only as a guideline. It is ultimately the responsibility of the laboratory to establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

PERFORMANCE CHARACTERISTICS

1. **Comparison:** A group of 40 sera ranging in BUN values from 5.9-108.5 mg/dL was assayed by the described BUN method and by a similar commercially available BUN reagent. Comparison of the results yielded a correlation coefficient of 0.994 and the regression equation was $y = 0.962x - 0.721$. (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T.)

2. **Precision:** Within-run precision was established by 20 assays on three different levels of commercial serum controls. Total precision values were obtained by assaying two commercial controls for five consecutive days.

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3. **Linearity:** 80 mg/dL. Performed according to NCCLS Guideline EP6-P.

REFERENCES

1. Sampson, E.J., Baird, M.A., Burtis, C.A., et al.: A coupled-enzyme equilibrium method for measuring urea in serum: Optimization and evaluation of the AACC Study Group on urea candidate reference method. Clin Chem, 26, 816-826, 1980.
2. Henry, R.J., Clinical Chemistry Principles and Techniques, 2nd Edition, Harper and Row, Hagerstown, NY, (1974), p. 511.
3. Tietz, N.W., Textbook of Clinical Chemistry, W.B. Saunders Company, Philadelphia, PA, (1986), p. 1270-1271.
4. Young, D.S., 3rd Ed., AACC Press, Washington DC (1990).

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INTENDED USE

Albumin reagent is used for the quantitative determination of albumin concentration in human serum.

INTRODUCTION

Albumin is the most abundant protein constituent of serum. It is synthesized in the liver and is noted for its ability of configuration changes. This steric affinity allows the albumin molecule to serve as a carrier of many substances such as bilirubin, fatty acids, uric acid, various drugs, and antibiotics. Albumin also functions in the maintenance of proper osmotic pressure.¹ Elevated serum albumin levels are associated with possible dehydration. Low serum albumin levels are indicative of potential malnutrition, liver disease, kidney disorders, and rheumatoid arthritis.² Earlier salt fractionation methods to determine serum albumin were too laborious to perform and have been replaced by azo dye methods. The use of bromocresol green in the reaction has become the preferred method because of its high specificity for albumin and its negligible interference for hemolysis, bilirubin, and salicylates.^{3,4}

PRINCIPLE

Serum albumin binds selectively to the dye bromocresol green at pH 4.2. The increase in absorbance of the resulting albumin-dye complex, read at 630 nm, is proportional to the albumin concentration.

REAGENT COMPOSITION

1. Bromocresol Green (BCG) 0.25 mM, buffer pH 4.0 ± 0.1, surfactant, non-reactive ingredients and stabilizers.
2. Standard: Bovine Albumin Fraction V with stabilizer (5 g/dl).

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
- CAUTION: *In vitro* diagnostic reagents may be hazardous. Handle in accordance with good laboratory procedures which dictate avoiding ingestion, and eye or skin contact. The reagent is an acid solution. Flush with water when contact occurs.
2. Specimens should be considered infectious and handled appropriately.
3. Avoid contact/ingestion. Reagent is an acid solution. Flush with water when contact occurs.

STORAGE AND STABILITY

1. Store the reagent at room temperature (15 - 30°C).
2. Albumin standard is stored refrigerated (2 - 8°C).

REAGENT DETERIORATION

The reagent should be a clear, yellow-green solution. If turbidity or precipitation has occurred the reagent should be discarded.

SPECIMEN COLLECTION

1. Test specimens should be serum and free from hemolysis.
2. Avoid excessive hemolysis since every 100 mg/dl of hemoglobin corresponds to about 100 mg/dl of albumin.
3. Albumin in serum is stable for one (1) week at room temperature (15 - 30°C) and approximately one (1) month when stored in the refrigerator.

INTERFERING SUBSTANCES

Ampicillin and other medications seriously interfere with the dye-binding properties of albumin (Young, et al.).⁵ It is also recommended that only standards and controls containing human albumin be employed with this procedure. The dye-binding properties of albumin from various species have been found to differ widely.

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1. Accurate pipetting devices
2. Test tubes/rack
3. Timer
4. Spectrophotometer

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1. Label test tubes Blank, Standard, Control, Patient, etc.
 2. Pipette 1.5 ml of reagent into each tube.
 3. Transfer 0.01 ml (10 µl) of sample to respective tubes and mix and allow to stand at room temperature for five (5) minutes.
 4. Zero spectrophotometer with the blank at 630 nm. (Wavelength range: 580 - 630 nm).
 5. Read and record absorbencies of all tubes.
- * ALTERNATIVE VOLUMES: 0.02 ml sample to 3.0 ml reagent.
* HTI MULTIPURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD

LIMITATIONS

1. The dye-binding properties of albumin, other than human, differ among species.
2. Samples with values above 8.0 g/dl should be diluted with 0.9% saline 1:1, re-run, and results multiplied by 2. Samples with results below 0.5 g/dl should be done electrophoretically.
3. Severely lipemic serums should have a serum blank.
 - a. Add 0.02 ml (20 µl) sample to 3.0 ml distilled water and read absorbance against distilled water at 630 nm.
 - b. Subtract the serum blank absorbance from the test absorbance and use the corrected absorbance in the calculations.

CALCULATIONS

Abs = Absorbance

Abs of unknown X Concentration of Standard = Albumin g/dl

Abs. of standard

Example:

If the absorbance of the unknown = 0.455, the absorbance of the standard is

0.705 and the standard concentration = 5.0 then:

0.455

0.705 x 5.0 = 3.23 g/dl

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established albumin values may be routinely used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction or procedural errors.

EXPECTED VALUES

3.5 - 5.3 g/dl

It is strongly recommended that each laboratory establish its range of expected values, since differences exist between instruments, laboratories, and local populations.

PERFORMANCE CHARACTERISTICS

1. Linearity: 0.5 - 8.0 g/dl

2. Sensitivity: Based on an instrument resolution of A = 0.001 the present method has a sensitivity of 0.005 g/dl.

3. Comparison: A comparison study performed between this method and another BCG method resulted in a correlation coefficient of 0.99 with a regression equation of $y = 0.96x + 0.1$.

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Intended Use

For the quantitative determination of ethyl alcohol in serum or whole blood. For *in vitro* diagnostic use only.


Clinical Significance

Alcohol analysis is most frequently requested for emergency room patients to aid in the differential diagnosis of central nervous system depression, which can be due to alcohol intoxication. It is also requested for forensic medicine purposes.

Method History

An enzymatic method for ethyl alcohol, using alcohol dehydrogenase and NAD, was first described by Bonnichsen et al.¹ Several modifications and improvements of this method have subsequently been reported. The present procedure is a modification of the original enzymatic procedure. Due to its speed and ease of use, enzymatic analysis of ethyl alcohol has become quite common in laboratories. But because of medicolegal issues, the analysis must also be highly accurate. Therefore, it is generally recommended that assays for ethyl alcohol be confirmed by another method such as gas chromatography.

Principle

Alcohol dehydrogenase (ADH) catalyzes the oxidation of ethanol to acetaldehyde with the concomitant reduction of NAD to NADH.² The change in absorbance at 340 nm is directly proportional to the alcohol concentration in the sample.

Reagents

Concentrations reagent to reconstituted reagent: Alcohol dehydrogenase (yeast) 200,000 U/L, NAD 2.6 mM, Buffer, pH 9.2 ± 0.1, surfactant, preservative.

Reagent Preparation

Reconstitute reagent with the volume of distilled or deionized water specified on the vial label. Swirl gently to dissolve.

Reagent Storage

1. Store reagent at 2-8°C.
2. Reconstituted reagent is stable for seven (7) days if tightly stoppered and stored at 2-8°C.

Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. The alcohol reagent is an IRRITANT. Handle with normal precautions. In case of contact rinse with plenty of water. Refer to MSDS for any hazard or safety information.
3. The reagent should not be used if moisture has penetrated the vial and caking has occurred or the reagent has an initial absorbance greater than 0.500 versus a water blank at 340 nm.
4. All specimens and controls should be handled in accordance with good laboratory practices using appropriate precautions as described in the CDC/NH Manual, "Biosafety in Microbiological and Biomedical Laboratories", 2nd Ed., 1988, HHS Publication No. (CDC) 88-8395.

Specimen Collection and Storage

1. The site of venipuncture should be disinfected only with aqueous disinfectants such as zephiran or merthiolate. Alcohol or other volatile disinfectants must not be used.
2. Whole blood may be collected using citrate, heparin or oxalate as the anticoagulant.³ However, fluoride used at 5 mg/dl of blood is the best preservative. Samples may be stored at 2-8°C for several days and may not show appreciable loss of alcohol, when well stoppered.⁴
3. Serum samples may be assayed without deproteinization.
4. Specimen collection should be carried out in accordance with NCCLS M29-T2.⁵ No method can offer complete assurance that human blood samples will not transmit

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Alcohol Reagent Set Алкоголь

Infection. Therefore, all blood samples should be considered potentially infectious.

Interferences

- The major interference in this assay is from alcohol used in cleansing the skin prior to venipuncture.
- A comparison of relative interference by various alcohols is listed below:

Alcohol reactive (after 5 minutes @ 30°C)	
Ethanol	100%
n-Butanol	38.5%
Isopropanol	6.6%
Methanol	0.0%
Ethylene Glycol	1.4%
Acetone	0.0%

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- Accurate pipetting devices
- Test tubes
- Timer
- Spectrophotometer able to read at 340 nm
- Water bath or Heating block (30°C or 37°C)
- Sodium Chloride, 0.9% (w/v). To prepare, dissolve 0.9 gm of reagent grade sodium chloride in distilled or deionized water and bring to a volume of 100 ml.
- Trichloroacetic Acid Solution (6.25% w/v) - if whole blood is to be measured.
- Centrifuge - if whole blood is to be measured.

Procedure (Automated-General)

Wavelength	340nm
Temperature	37°C or 30°C
Assay Type	Endpoint
Sample/Reagent Ratio	1:201
Reaction Direction	Increasing
Incubation Time	5 minutes
Low Normal	0 mg/dl
High Normal	100 mg/dl

Application parameters for various automated instruments are available. Please contact the manufacturer's Technical Service Department for specific information.

Procedure (Manual)

- Reconstitute reagent according to instructions.
- Pipette 1.0ml of reagent into the required test tubes.
- Bring the reagent in the tubes to the selected temperature (30°C is recommended for manual assay) and add 5 µl of the standard, test and control samples to the respective tubes. To the blank tube, add 5 µl of 0.9 sodium chloride solution.
- Mix gently. Incubate for 5 minutes.
- At the end of the incubation period, read and record the absorbance of all tubes against the reagent blank at 340 nm. Use this absorbance reading in the calculations.

NOTE: To measure whole blood samples the specimen must be deproteinized as follows:

- Pipette 1.8ml Trichloroacetic Acid Solution into a centrifuge tube.
- While swirling tube, slowly add 0.2ml sample.
- Stopper the tube and mix vigorously. Allow to stand at room temperature for approximately five minutes.
- Centrifuge (2000 rpm) for five minutes to obtain a clear supernatant.
- This supernatant may be assayed as above. However, to account for the dilution of the specimen it is recommended to use a 1:10 dilution of standard for calibration of the assay.

Limitations

Alcohol concentrations greater than 400 mg/dl should be diluted 1:1 with normal saline.

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Multiply results by 2 to compensate for dilution.

Calibration

Use an NIST-traceable Ethanol Standard or serum calibrator. The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the procedure should be re-calibrated.

Calculation

Ethanol conc. (mg/dl) = $\frac{\text{Abs. Sample}}{\text{Abs. Std.}} \times \text{conc. of Std. Abs. Std.}$

Example Calculation: Abs of patient = 0.450, Abs. of standard = 0.500, Conc. of Standard = 100

$$(0.450 / 0.500) \times 100 = 90 \text{ mg/dl}$$

Quality Control

Controls are recommended to monitor the performance of the assay, providing a constant screening of the instrument, reagents and techniques. Commercially available control material with established values for alcohol may be used. The controls should be run at least with every working shift in which ethanol assays are performed. It is recommended that each laboratory establish their own frequency of control determination.

Expected Values

Patients abstaining from alcohol have non-detectable levels by most enzymatic or chromatographic methods.⁸ Below are the blood alcohol ranges shown with the corresponding percentage of subjects designated as intoxicated. In a study of over 6000 subjects, indications are very few (about 4%) people are intoxicated at blood levels of 0.05%.⁹

Blood Alcohol Levels	Percentage Diagnosed as Intoxicated
0.1 - 0.15%	63%
0.15 - 0.2%	89%
0.2 - 0.25%	95%

Above based on the averaged results of several studies.⁸ concentrations are usually about 16% higher than blood samples from which they are derived.⁸

Performance

- Assay Range: 0 - 400 mg/dl
- Comparison: Testing performed on 100 specimens compared with a similar method for ethanol yielded a correlation coefficient of 0.999 with a regression equation of $y = 1.01x - 1.01$. Sample values ranged from 2 to 412. ($Sy, x = 12.84$)
- Precision: Precision studies were performed following a modification of the guidelines contained in NCCLS document EP5-T2.¹¹

Within Run			Run to Run		
Mean	S.D.	C.V.%	Mean	S.D.	C.V.%
40.3	0.72	1.78%	40.6	1.36	3.34%
93.3	1.04	1.11%	89.0	3.27	3.67%
219.8	2.17	0.99%	224.7	3.71	1.65%

- Sensitivity: The sensitivity for the Alcohol reagent was investigated by reading the change in absorbance at 340 nm for a saline sample, and serums with known concentrations. Ten replicates of each sample were performed. The results of this investigation indicated that on the analyzer used, the Alcohol reagent showed little or no reagent drift on a zero sample. Also, that an absorbance change of 0.006 was approximately equivalent to one mg/dl of ethanol.

References

1.2.

- Bonishen, R.K., Theorell, H., Scand. J. Clin. Lab. Invest. 3:58, 1951.
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues, In Methods of Biochemical Analysis, Vol. VII, D. Gluck, Editor, Interscience, New York, pp. 217-251, 1957.
- Jones, D., Gerber, L.P., Drell, W., Clin. Chem. 16:402, 1970.
- Hepler, O.E., Manual of Clinical Laboratory Methods, 4th Ed., C.C. Thomas, Springfield, IL, p. 329, 1949.
- NCCLS document "Protection of Laboratory Workers from Infectious Disease

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Transmitted by Blood, Body Fluids, and Tissue", 2nd Ed., 1991.
Data obtained by the manufacturer.
Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACCC Press,
Washington, D.C.; Supplement No., 1991.
Dubowski, K.M., Laboratory Management 4.27, 1982.
American Medical Association: Alcohol and the Impaired Driver, 1966.
Shoemaker, M.J., Pathologist, 26, 1985.
NCCLS document "Evaluation of Precision Performance of Clinical
Chemistry Devices", 2nd Ed., 1992.

Manufactured by:
USA, HTI Medical 109 production Rd Walpole, MA 02081

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INTENDED USE

For the quantitative determination of alanine aminotransferase in serum used in routine examination and monitoring of therapy and relapses.


INTRODUCTION

The enzyme alanine aminotransferase is widely reported in a variety of tissue sources. The major source of ALT is of hepatic origin and has led to the application of ALT determinations to the study of hepatic diseases. Elevated serum levels are found in hepatitis, cirrhosis, and obstructive jaundice. Levels of ALT are only slightly elevated in patients following a myocardial infarction.

UV methods for ALT determination were first developed by Wroblewski and LaDue in 1956. The method was based on the oxidation of NADH by lactate dehydrogenase (LDH). In 1980, the International Federation of Clinical Chemistry recommended a reference procedure for the measurements of ALT based on the Wroblewski and LaDue procedure. The ALT reagent conforms to the formulation recommended by the IFCC.

PRINCIPLE

The enzymatic reaction sequence employed in the assay of ALT is as follows:

The pyruvate formed in the first reaction is reduced to lactate in the presence of lactate dehydrogenase and NADH. The activity of ALT is determined by measuring the rate of oxidation of NADH at 340 nm. Endogenous sample pyruvate is converted to lactate by LDH during the lag phase prior to measurement.

REAGENTS

ALT Liquid Reagents 1 and 2 come in separate containers, and both reagents are clear, colorless liquid in ready to use format. After combining ALT Liquid R1 (Buffer Reagent) and ALT Liquid R2 (Co-Enzyme) the working reagent contains:

L-Alanine	500 mmol/L
LDH	>1200 U/L
Tris Buffer, pH 7.5	100 mmol/L
2-Oxoglutarate	15 mmol/L
NADH (Disodium salt)	0.18 mmol/L

Stabilizers and Preservatives

WARNINGS AND PRECAUTIONS

Normal precautions exercised in handling laboratory reagents should be followed. The reagents contain sodium azide, which may be toxic if ingested. Sodium azide may also react with lead and copper plumbing to form highly explosive metal azides. Refer to Material Safety Data Sheet for any updated risk, hazard, or safety information.

REAGENT PREPARATION

The working reagent is prepared by mixing five (5) volumes of R1 with one (1) volume of R2 in a disposable container. Example: 25 ml R1 + 5 ml R2

REAGENT STORAGE

Reagents are stable until the expiration date on their respective labels, when properly stored at 2-8°C and protected from light. Reagents should appear clear and colorless.

REAGENT DETERIORATION

- Discard if either appears cloudy or contains particulate matter.
- The working reagent is stable for 2 weeks at 2-8°C. The working reagent should be discarded if the initial absorbance, read against distilled water at 340 nm, is below 1.000.

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- Spectrophotometer capable of absorbance reading at 340 nm and 1 cm light path
- Constant temperature block or bath, 37°C, or temperature controlled cuvette well
- Accurate pipetting devices
- Test tubes
- Interval timer

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SPECIMEN COLLECTION AND STORAGE

Non-hemolyzed serum is the specimen of choice. Whenever possible specimens should be separated and analyzed on the day of collection. Store serum in stoppered tubes. About 10% ALT is lost 3 days at 4°C and in 1 day at 25°C.

INTERFERING SUBSTANCES

Hemolysis must be avoided as the concentration of ALT in red cells is roughly 5 times that of serum. Bilirubin levels up to 40 mg/dl, and triglyceride levels up to 2000 mg/dl, show no interference in this test. Certain drugs and other substances are also known to affect ALT values.

MANUAL PROCEDURE

1. Prepare ALT working reagent according to instructions.
2. Pipette 1.0 mL of working reagent into tubes labeled "controls", "patient(s)", etc.
3. Pre-incubate all tubes at 37°C for at least five minutes.
4. Zero spectrophotometer at 340 nm with distilled water.
5. Add 100 µL (0.10 mL) serum to its respective tube, mix gently and turn to a thermo cuvette.
6. Read and record absorbance at 1 minute. Continue incubating at 37°C and record absorbance again at 2 and 3 minutes. Rate should be constant.
7. Determine the average absorbance per minute (AA/min), multiply by factor 1768 for results in U/L.
8. Repeat the procedure for each sample.

NOTE: If cuvette is not temperature controlled, incubate samples at 37°C between readings.

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QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established ALT values may be routinely used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedure errors.

CALIBRATION

ALT activity is based on the "micromolar extinction coefficient" of NADH at 340 nm (see "Results" section). The instrument manufacturer's calibration guidelines should be followed to calibrate your analyzer. Assaying the ALT contents of a control serum with known ALT values can be used to assure instrument calibration has been performed correctly.

RESULTS

Values are derived based on the "absorptivity micromolar extinction coefficient" of NADH at 340 nm (0.00622). Units per liter (U/L) of ALT/GPT activity is that amount of enzyme, which oxidizes one [µmol] of NADH per minute.

$$U/L = \frac{AA/min}{Absorptivity} \times \frac{Total Volume}{Sample Volume}$$

$$U/L = \frac{AA/min}{0.00622} \times \frac{1.100}{0.100}$$

$$U/L: AA/min \times 1768$$

LIMITATIONS

If the AA/min is greater than 0.342, dilute 1 part sample with 9 parts isotonic saline and re-assay. Multiply the result by 10. ALT values for neonatal patients have not been established with this procedure. Grossly icteric or turbid specimen may require the use of a sample blank.

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PERFORMANCE CHARACTERISTICS

1. Comparison: A group of 128 sera ranging in ALT activity from 7 - 625 U/L was assayed by the described ALT method and by a similar commercially available ALT reagent. Comparison of the results yielded a correlation coefficient of 0.999 and the regression equation was $y = 0.960x + 3.2$. (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T.)

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3. Linearity: Linear to 500 U/L at 37°C. Performed according to NCCLS Guideline EP6-P.

4. Sensitivity: Based on an instrument resolution of $A = 0.001$, the method presented shows a sensitivity of 1.8 U/L.

REFERENCES

1. Henry, J.B. *Clinical Diagnosis and Management by Laboratory Methods*, W.B. Saunders and Co., Philadelphia, PA, p-332-335 (1974).
 2. Wroblewski, F. and LaDue, J.S. *Proc. Soc. Exper. Biol. and Med.* 91:569 (1956).
 3. International Federation of Clinical Chemistry, *J. Clin. Chem. Clin. Bio.* 18:5231 (1980).
 4. Bergmeyer, H.U. *Principles of Enzymatic Analysis*. Verlag Chemie, 1978.
 5. Young D.S. *Effects of drugs on clinical laboratory tests*. AACCPress, Washington D.C., 1990.
 6. Henry J.B. *Clinical Diagnosis and Management by Laboratory Methods*, 17th ed. WB Saunders Co., 1964, p 1437.
- AS24: 10/01
AA/Min Total Volume
Absorbivity Sample Volume
AA/Min

Manufactured by: HTI Medical
109 Production Road
Walpole, MA 02081
USA

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Alpha-1 Microglobulin

Type: Reagent
Method: Immunoturbidimetric
Format: Liquid
Linearity: 1.0 - 137.0 mg/L (or value of highest calibration point)
Urine: 0.2 - 34.0 mg/L (or value of highest calibration point)
Sample: Serum, Plasma, Urine

KIT COMPOSITION:

Note: Typically manufactured only by request.

R1: 1 X 50 mL buffer
R2: 1 X 50 mL antibody
Tests / Kit:
277 serum/plasma (Hitachi 717)
400 urine (Hitachi 717)
(Can Vary Based On Analyzer Model)

INTENDED USE:

For the quantitative determination of human Alpha-1 Microglobulin in human urine, serum, and plasma by immunoturbidimetric assay.
FOR IN VITRO DIAGNOSTIC USE.

Principle

Endpoint Measuring reaction antigen-antibody.

Reagent Composition

Latex Solution
Glycin Buffer (pH 8.2)
Latex particles sensitized with goat antibodies to AMI (0.25%)
Sodium azide (0.95 g / l)
Buffer
NaCl (9 g / l)
Detergent (0.01%)
Sodium azide (0.95 g / l)

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STORAGE AND HANDLING:

All reagents should be stored refrigerated (2-8°C). Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for one year from the date of manufacture, as indicated by the expiration date on the package and bottle labels. Opened reagents can be used for one month if stored at 2-8°C. It is recommended that R-2 be gently inverted once a week after it is opened.

PERFORMANCE:

Precision - Within Run

Serum:	Sample II	Sample III
Sample I	N = 10	N = 10
N = 10	Mean = 32.1	Mean = 66.1
Mean = 7.1	S.D. = 0.263	S.D. = 0.540
S.D. = 0.145	C.V. = 0.82%	C.V. = 0.82%
C.V. = 2.04%		

Urine:

Sample I	Sample II	Sample III
N = 10	N = 10	N = 10
Mean = 0.7	Mean = 3.3	Mean = 17.4
S.D. = 0.029	S.D. = 0.012	S.D. = 0.059
C.V. = 3.93%	C.V. = 0.37%	C.V. = 0.34%

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ALPHA-FETOPROTEIN (AFP) ENZYME IMMUNOASSAY TEST KIT

Enzyme Immunoassay for the Quantitative Determination of Alpha-Fetoprotein (AFP) in Human Serum

Store at 2 to 8°C.

INTENDED USE

For the quantitative determination of the Cancer Antigen AFP concentration in human serum.

INTRODUCTION

Alpha-fetoprotein (AFP) is a glycoprotein with a molecular weight of approximately 70,000 daltons. AFP is normally produced during fetal and neonatal development by the liver, yolk sac, and in small concentrations by the gastrointestinal tract. After birth, serum AFP concentrations decrease rapidly, and by the second year of life and thereafter only trace amounts are normally detected in serum.

Elevation of serum AFP to abnormally high values occurs in several malignant diseases, most notably nonseminomatous testicular cancer and primary hepatocellular carcinoma. In the case of nonseminomatous testicular cancer, a direct relationship has been observed between the incidence of elevated AFP levels and the stage of disease. Elevated AFP levels have also been observed in patients diagnosed with seminoma with nonseminomatous elements, but not in patients with pure seminoma. In addition, elevated serum AFP concentrations have been measured in patients with other noncancerous diseases, including ataxia telangiectasia, hereditary tyrosinemia, neonatal hyperbilirubinemia, acute viral hepatitis, chronic active hepatitis, and cirrhosis. Elevated serum AFP concentrations are also observed in pregnant women. Therefore, AFP measurements are not recommended for use as a screening procedure to detect the presence of cancer in the general population.

Principle of the test

The AFP ELISA test is based on the principle of a solid phase enzyme-linked immunosorbent assay. The assay system utilizes a goat anti-AFP antibody directed against intact AFP for solid phase immobilization (on the microtiter wells). A monoclonal anti-AFP antibody conjugated to horseradish peroxidase (HRP) is in the antibody-enzyme conjugate solution. The test sample is allowed to react first with the immobilized rabbit antibody for 30 minutes. The wells are washed to remove any unbound antigen. The monoclonal-HRP conjugate is then reacted with the immobilized antigen for 30 minutes at room temperature resulting in the AFP molecules being sandwiched between the solid phase and enzyme-linked antibodies. The wells are washed with water to remove unbound labeled antibodies. A solution of TMB Reagent is added and incubated for 20 minutes, resulting in the development of a blue color. The color development is stopped with the addition of Stop Solution changing the color to yellow. The concentration of AFP is directly proportional to the color intensity of the test sample. Absorbance is measured spectrophotometrically at 450 nm.

Reagents

Materials provided with the kits:

- Goat anti-AFP coated microtiter plate with 96 wells.
- Zero Buffer, 13 ml.
- Reference standard set, contains 0.5, 20, 50, 150, and 300 ng/ml (WHO, 72/225) AFP, lyophilized.
- Enzyme Conjugate Reagent, 18 ml.
- TMB Reagent (One-Step), 11 ml.
- Stop Solution (1N HCl), 11 ml.

Materials required but not provided:

- Precision pipettes: 20 µl, 100 µl, and 150 µl and 1.0 ml.
- Disposable pipette tips.
- Distilled water.
- Vortex mixer or equivalent.
- Absorbent paper or paper towel.
- Graph paper.
- Microtiter plate reader.

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Serum should be prepared from a whole blood specimen obtained by acceptable medical techniques. This kit is for use with serum samples without additives only.

Storage of Test Kit and Instrumentation:

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Unopened test kits should be stored at 2-8°C upon receipt and the microtiter plate should be kept in a sealed bag with desiccants to minimize exposure to damp air. Opened test kits will remain stable until the expiration date shown, provided it is stored as described above. A microtiter plate reader with a bandwidth of 10 nm or less and an optical density range of 0-2 OD or greater at 450 nm wavelength is acceptable for use in absorbance measurement.

Reagent Preparation:

1. All reagents should be brought to room temperature (18-25°C) before use.
2. Reconstitute each lyophilized standard with 1.0 ml distilled water. Allow the reconstituted material to stand for at least 20 minutes and mix gently.
3. Reconstituted standards will be stable for up to 30 days when stored sealed at 2-8°C.

ASSAY PROCEDURE

1. Secure the desired number of coated wells in the holder.
2. Dispense 20 µl of standard, specimens, and controls into appropriate wells.
3. Dispense 100 µl of Zero Buffer into each well.
4. Thoroughly mix for 30 seconds. It is very important to have a complete mixing in this setup.
5. Incubate at room temperature (18-25°C) for 30 minutes.
6. Remove the incubation mixture by flicking plate content into a waste container.
7. Rinse and flick the microtiter wells 5 times with distilled or deionized water.
8. Strike the wells sharply onto absorbent paper or paper towels to remove all residual water droplets.
9. Dispense 150 µl of Enzyme Conjugate Reagent into each well. Gently mix for 10 seconds.
10. Incubate at room temperature for 30 minutes.
11. Remove the incubation mixture by flicking plate contents into a waste container.
12. Rinse and flick the microtiter wells 5 times with distilled or deionized water. (Please do not use tap water.)
13. Strike the wells sharply onto absorbent paper to remove residual water droplets.
14. Dispense 100 µl TMB Reagent into each well. Gently mix for 10 seconds.
15. Incubate at room temperature for 20 minutes.
16. Stop the reaction by adding 100 µl of Stop Solution to each well.
17. Gently mix for 30 seconds. **It is important to make sure that all the blue color changes to yellow color completely.**
18. Read optical density at 450 nm with a microtiter reader within 15 minutes.

CALCULATION OF RESULTS

1. Calculate the average absorbance values (A450) for each set of reference standards, control, and samples.
2. Construct a standard curve by plotting the mean absorbance obtained for each reference standard against its concentration in ng/ml on linear graph paper, with absorbance on the vertical (y) axis and concentration on the horizontal (x) axis.
3. Using the mean absorbance value for each sample, determine the corresponding concentration of AFP in ng/ml from the standard curve.

EXAMPLE OF STANDARD CURVE

Results of a typical standard run with optical density readings at 450 nm shown in the Y axis against AFP concentrations shown in the X axis. This standard curve is for the purpose of illustration only, and should not be used to calculate unknowns. Each user should obtain his or her own data and standard curve.

AFP (ng/ml)	Absorbance (450 nm)
0	0.012
5	0.127
20	0.465
50	0.952
150	2.150
300	2.932

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
EXPECTED VALUES AND SENSITIVITY

In high-risk patients, AFP values between 100 and 350 ng/ml suggest a diagnosis of hepatocellular carcinoma; and levels over 350 ng/ml usually indicate the disease. Approximately 97% of the healthy subjects have AFP levels less than 8.5 ng/ml. It is recommended that each laboratory establish its own normal range. The minimum detectable concentration of AFP by this assay is estimated to be 2.0 ng/ml.


Limitations of the Procedure

Reliable and reproducible results will be obtained when the assay procedure is carried out with a complete understanding of the package insert instructions and with adherence to good laboratory practice.

The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.

Serum samples demonstrating gross lipemia, gross hemolysis, or turbidity should not be used with this test.

The results obtained from the use of this kit should be used only as an adjunct to other diagnostic procedures and information available to the physician.

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Intended Use

This reagent is intended for the enzymatic measurement of ammonia in plasma. For *in vitro* diagnostic use only.


Clinical Significance

The bulk of ammonia in the body is generated in the gastrointestinal system by action of bacterial enzymes on the contents of the colon and from hydrolysis of glutamine. It is removed in the liver and converted to urea through a series of enzymatic reactions in the Krebs-Henseleit cycle. Among other conditions, advanced liver disease and hepatic encephalopathy result in elevated levels of ammonia in blood. Hyperammonemia is also common in inherited deficiencies of the enzymes involved in the conversion of ammonia to urea.

The determination of ammonia is very useful in the diagnosis and prognosis of Reye's Syndrome. Elevated blood ammonia exerts toxic effects on the central nervous system [1,2].

Test Summary

The enzymatic determination of ammonia allows a direct measurement of the compound in the plasma which avoids the long and laborious methods of separation employed in older methodologies. The enzymatic assay gives a highly sensitive and specific method. The assay is based on the following reaction [2,3]:

Ammonia reacts with α -ketoglutarate (α -KG) and reduced nicotinamide adenine dinucleotide phosphate (NADPH) to form L-glutamate and NADP in a reaction catalyzed by glutamate dehydrogenase (GLDH) (L-glutamate: NAD(P)⁺ oxidoreductase (deaminating), EC 1.4.1.3). The amount of NADPH oxidized is, on a molar basis, equal to the content of ammonia in the sample. The reaction can be followed by the decrease in absorbance at 340nm. The reagent is provided in two separate vials. Keeping the components of the reagent separated until time of assay increases their stability after reconstitution. The use of NADPH in place of NADH minimizes interference by such components of plasma as pyruvate and lactate dehydrogenase.

Reagents

Ammonia Substrate/Reagent 1
Buffer 100 mmol/L
EDTA 2 mmol/L
 α -Ketoglutarate 3.4 mmol/L
Adenosine diphosphate 0.5 mmol/L
NADPH 0.3 mmol/L
Fillers and Stabilizers
pH=8.6 ± 0.1
Ammonia Enzyme/Reagent 2
Buffer 100 mmol/L
EDTA 2 mmol/L
Adenosine diphosphate 0.5 mmol/L
GLDH 400 KU/L
Fillers and Stabilizers
pH=7.8 ± 0.1

Reagent Preparation

Use only ammonia-free water, freshly deionized or distilled. Dissolve Reagent 1 and Reagent 2 with the volume of water specified on the vial label. Maintain solutions tightly capped.

Reagent Storage and Stability

The dry reagents in the unopened vials are stable until the expiration date on the vial labels. The reconstituted reagents are stable for at least 15 days stored at 2°C to 8°C.

Precautions

1. Avoid ammonia contamination from the air, water, and glassware. Ammonia contamination can be checked by assaying the water used with these reagents. Run a blank assay by substituting the water used for the sample.

2. Through-traffic and smoking must be avoided in the patient's room and in the laboratory where the assay is performed. The phlebotomist should be a

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nonsmoker.

If the patient is a smoker, wash site of venipuncture. Blood should be drawn in a room where no smoking is permitted.

3. Do not use reagent if the absorbance of Reagent 1 read at 340 nm against a blank of water is less than 1,200.

4. Avoid contaminating reagents. If reagent shows microbial contamination, as indicated by turbidity, do not use.

5. WARNING: This product contains < 0.1% sodium azide. Sodium azide may react with lead or copper plumbing to form explosive compounds. When disposing of this product through plumbing fixtures, flush with large amounts of water to prevent azide build up.

Specimen Collection and Storage

EDTA plasma is the specimen of choice. The use of heparin as an anticoagulant is not recommended. Collect blood from a stasis-free vein into an EDTA evacuated tube, release residual vacuum in the tube; mix gently, place on ice and deliver to the laboratory without delay. Separate the plasma from the cells immediately. Do not use hemolyzed samples. The analysis should be performed within 30 minutes. A maximum of 2 hours delay with the plasma on ice is permissible.

Interferences

The major interference for this assay is from contamination by ammonia in the air and water. Analytical and physiological variables including drugs and other substances which influence ammonia concentrations have been listed by Young 4, 5, 6.

Materials Provided

Ammonia Reagents in dry powder form in vials to be reconstituted as described above.

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1. Accurate pipetting devices
2. Matched cuvettes, square, 1 cm light path.
3. Ammonia standards or calibrators
4. Spectrophotometer able to read at 340 nm.
5. Constant temperature device is not needed. However, use the same temperature for assay of controls and samples.

Procedure

Wavelength: 340nm
Temperature: 25°C
Cuvettes (square): 1 cm light path.
Blank: Water
Reagent 1: 1 mL
Sample: 0.2 mL

1. Mix. Incubate for 4 minutes. Read absorbance at 340nm with instrument set to zero absorbance with the water blank. This is reading R1, which must be corrected to compensate for the volume addition of Reagent 2.

$R1 \times 0.96 = R1c$

R1c is used to calculate the ΔA below.

2. Add 0.05ml Reagent 2. Mix. Incubate. After 5 minutes read absorbance again. This is reading R2. Calculate the change in absorbance, ΔA . Use this ΔA value in the calculation below.

$\Delta A = R1c - R2$

Note: Run a blank assay using water in the place of the sample to check for contamination. The value obtained for the blank should then be subtracted from the value found for the sample.

Limitations

Samples with ammonia concentrations exceeding 600 $\mu\text{mol/L}$ ($\Delta A > 0.600$) should be assayed again after dilution with an equal volume of distilled or deionized water. Multiply result by 2.

Calibration

Assay values can be obtained with this test procedure by using a factor or by using ammonia standards. (Refer to the Calculations section below.)

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Ammonia Reagent Set Аммиак

Commercially available ammonia standards or calibrators may be used. Ammonium sulfate standards also can be obtained from the College of American Pathologists or can be prepared from desiccated reagent grade ammonium sulfate.

Calculations

1 Factor
 $\Delta A \times 1.25 \times 1000 = \Delta A \times 1005 = \text{umol/L ammonia in sample}$
 $6.22 \times 1 \times 0.2$
 Where: 1.25 = Total volume of assay (mL)
 1000 = Conversion to liter volume
 6.22 = Millimolar extinction coefficient of NADPH at 340nm
 1 = Light path (cm)
 0.2 = Volume of the sample (mL)
 The factor 1005 multiplied by ΔA equals the ammonia concentration in the sample in umol/L.
 2. Standard
 Concentration of Standard $\times \Delta A$ Sample = umol/L ammonia
 ΔA of Standard in Sample

Sample Calculations

1 Factor
 If ΔA of sample is 0.050 then, $0.050 \times 1005 = 50 \text{ umol/L ammonia in sample}$
 2. Standard
 Concentration of ammonia standard = 500 umol/L
 ΔA of ammonia standard = 0.496
 ΔA of sample = 0.115
 500
 $\frac{500}{0.496} \times 0.115 = 116 \text{ umol/L ammonia in sample}$
 0.496

Quality Control

Controls are recommended to monitor the performance of the assay, providing a constant screening of the instrument, reagents and technique.

Expected Values

The expected range reported for the enzymatic procedure described is 11 umol/L to 35 umol/L. It is recommended that each laboratory establish its own reference range.

Performance

1. Linearity: 600 umol/L.
 2. Comparison: Sixty patient samples were divided into two respective aliquots and stored in an ice bath. Analyses were performed at the same time in two separate locations using this reagent on the Cobas FaraTM and the reagent by DuPont on the DuPont ACATM. The samples ranged in ammonia concentrations from 8 umol/L to 347 umol/L. The correlation coefficient was 0.999 and the regression equation was: $y = 0.991x - .561$.

Precision

Within Run n=10

Mean	S.D.	C.V.%
94	1.3	1.37
209	1.47	0.70
411	0.92	0.22

References

1. Textbook of Clinical Chemistry Edited by N.W. Tietz, W.B. Saunders Company, Philadelphia, p. 1409, 1986.
2. Ratcliff, C.R. and Hall, F.F. in Selected Methods of Clinical Chemistry, volume 9, p. 85. Edited by Willard R. Faulkner and Samuel Meites. American Association for Clinical Chemistry. Washington, D.C., 1982.
3. Da Fonseca-Wollheim F., J. Clin. Chem. Clin. Biochem. 11, 421, 1973.
4. Young, D.S., Effects of Pre-analytical Variables on Clinical Laboratory Tests, First Edition, AACC Press, Washington, D.C., 3.20-3.21, 1993.
5. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Third Edition, AACC Press, Washington, D.C., 3.30-3.32, 1990.
6. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, 1991 Supplement to the Third Edition, AACC Press, Washington, D.C., 3.9-3.10, 1991.

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
INTENDED USE

For the quantitative determination of α -amylase activity in serum using a manual or automated procedure.

PRINCIPLE OF THE METHOD

α -Amylase hydrolyzes the 2-chloro-4-nitrophenyl- α -D-maltotriose (CNP3) to release 2-chloro-4-nitrophenol (CNP) and form 2-chloro-4-nitrophenyl- α -D-maltoside (CNP2), maltotriose (G3) y glucose (G) according to the following reaction

The rate of 2-chloro-4-nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of α -amylase present in the sample

CLINICAL SIGNIFICANCE

α -Amylase (AMS) is an enzyme that helps to digest the glycogen and the starch. It is produced mainly by exocrine pancreas and salivary glands. This determination is made mainly in diagnosis or to control diseases of the pancreas as acute or chronic pancreatitis. It can also reflect biliary or gastrointestinal disease and other upheavals^{2,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data

REAGENT PREPARATION

No preparation is required

REAGENT COMPOSITION

MES 100 mmol/L, pH 6.0
CNP3 2.25 mmol/L
Sodium chloride 350 mmol/L
Calcium acetate 6 mmol/L
Potassium thiocyanate 900 mmol/L
Sodium azide 0.95 g/L

REAGENT STORAGE AND STABILITY

HTI α -Amylase Reagent stored unopened at 2°C to 8°C is stable until the expiration date showed on the bottle label. Once opened, HTI α -Amylase Reagent is stable for 30 days, unless the expiration date is exceeded.

DO NOT FREEZE

SPECIMEN COLLECTION AND HANDLING

1. The test can be performed on serum, plasma. For serum, blood is drawn into a tube which does not contain anticoagulant and allow clotting. The serum is then separated from the clot. A maximum limit of two hours from the time of collection is recommended.
2. Separated serum or plasma should not remain at room temperature longer than 8 hours. If assays are not completed within 8 hours, serum and plasma should be stored at 2°C to 8°C. If assays are not completed within 48 hours, or the separated sample is to be stored beyond 48 hours, samples should be frozen at -15°C to -20°C. Frozen samples should be thawed only once. Analyte deterioration may occur in samples that are repeatedly frozen and thawed.
3. It is recommended that urine assays be performed within 2 hours of collection. Due to the instability of amylase in acidic urine, the pH of the specimen should be adjusted to the alkaline range and stored at 4°C. For timed specimens, the collection container should be kept in the refrigerator or on ice during the timed period.
4. For plasma, add whole blood directly into a tube containing anticoagulant. Acceptable anticoagulants are listed in the "LIMITATIONS" section.

CALIBRATION

Calibration is not required.

MATERIALS NEEDED BUT NOT SUPPLIED WITH REAGENT KIT

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Citrate and EDTA were found to be incompatible with this method.
2. The anticoagulants Ammonium Heparin, Sodium Heparin and Lithium Heparin were found to be compatible with this method.

INTERFERENCE

1. Samples showing evidence of hemolysis should not be used.
2. Lipemic samples >3+ should be ultra-centrifuged and the analysis performed on the infranate.
3. Pyruvate at a level of 2 mg/dL may cause decreased results.
4. On this method, refer to the work of Young for a review of drug and comprehensive list of substances effect on α-Amylase level.

EXPECTED VALUE

Serum : 25-125 IU/L.

PRECAUTIONS

1. For in vitro diagnostic use only.
2. Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
3. Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.

MANUAL PROCEDURE

1. Prepare the working reagent according to instructions.
2. Pipette 1.0 ml of reagent into tubes labeled "Control", "Patient" etc.
3. Pre-warm all tubes at 37°C for at least 3 minutes.
4. Zero spectrophotometer with water at 405 nm. Add 0.025ml (25 µl) of sample and read after 60 seconds.
5. Continue reading every 60 seconds for 2 minutes.
6. Determine the mean absorbance difference per minute (ΔAbs./min.).
7. Multiply the ΔAbs./min. by 3178 to obtain result in IU/L.

8. TC MATRIX PROCEDURE

TEST NAME: AMYLASE-CNPG3	R1:	200
TEST NO.	R2:	0
FULL NAME: α-amylase	SAMPLE VOLUME:	5
REFERENCE NO.	R1 BLANK:	/
ANALY TYPE: Kinetic	MIX REAG. BLANK:	/
PRI. WAVE: 405nm	CONCENTRATION:	/
SECON. WAVE	LINEARITY LIMIT:	/
TREND: Increase	SUBSTRATE LIMIT:	/
REACT. TIME: 4 - 12	FACTOR: 3178	
INCUBATE TIME: /	PROZONE CHECK:	/

9.	UNIT: IU/L	Q1: / Q2: / Q3: / Q4: /
	PRECISION: Integer	PC: / ABS: /

CALIBRATION

The procedure is standardized by means of the millimolar absorptivity of 2-chloro-p-nitrophenol that is 12.9 at 405 nm under the test conditions described.

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1000 = Conversion of U/ml to U/L
MMA = Millimolar absorptivity of
2-chloro-p-nitrophenol (12.9)
SV = Sample volume (0.025 ml)
LP = Light path (1 cm)

$$\Delta\text{Abs./min} \times 1.025 \times 1000 = \Delta\text{Abs./min} \times 3178 = \text{U/L } \alpha\text{-amylase}$$

$$12.9 \times 0.025 \times 1.0$$

Example: If $\Delta\text{Abs./min} = 0.03$, then $0.03 \times 3178 = 95 \text{ U/L}$
NOTE: To convert to SI Units (kU/L) multiply the U/L value by 16.67.

PERFORMANCE CHARACTERISTICS

Analytical Range: 13-2400 U/L
For α -amylase analyte by HTI α -Amylase Reagent, this method has been demonstrated to be linear from 13-2400 U/L.

Accuracy: Comparison study was performed on Beckman CX System from 40 samples. Beckman Coulter α -Amylase Reagent was used to compare with HTI α -Amylase Reagent. The results of this study in yield a correlation coefficient of 0.999 with a regression equation of $y=1.0253x-7.0454$.

Precision: Within Run precision for HTI α -Amylase Reagent Set was determined following a modification of NCCLS EP5-A. Two commercial human serum were assayed on TC Matrix for 25 times.

Sample	Sample 1	Sample 2
N	25	25
Mean (mg/dl)	61	376
Standard Deviation	0.97	5.43
Coefficient of Variation	1.6	1.4

Run-Day precision for HTI α -Amylase Reagent was determined following a modification of NCCLS EP5-A. Two commercial human serum were assayed on TC Matrix five times per day for five days for the total of 25 values.

Sample	Sample 1	Sample 2
N	25	25
Mean (mg/dl)	62	391
Standard Deviation	1.13	7.31
Coefficient of Variation	1.8	1.9

REFERENCES:

1. Richter R. 1969. Clinical chemistry. Theory and Practice, A.P. New York.
2. Caraway W.T. 1959. Am. Clin. 32:97.
3. Tietz, N.W., "Specimen Collection and Processing: Sources of Biological Variation," Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders, Philadelphia, PA (1994)
4. National Committee for Clinical Laboratory Standards. Approved Guideline, NCCLS publication C28-A, Villanova, PA (1994).
5. Henry, J. B., ed., Clinical Diagnostics and Management by Laboratory Methods, 18th Edition, W.B. Saunders, Philadelphia.
6. Tietz, N.W., ed., Clinical Guide to Laboratory Tests, 2nd Edition, W.B. Saunders, Philadelphia, PA (1990)
7. National Committee for Clinical Laboratory Standards, Method Comparison and Bias Estimation Using Patient Samples, Tentative Guideline, NCCLS Publication EP9-T, Villanova, PA (1993)
8. National Committee for Clinical Laboratory Standards, Precision Performance of Clinical Chemistry Devices, Tentative Guideline, 2nd Edition, NCCLS publication EP5-T2, Villanova, PA (1992)
9. National Committee for Clinical Laboratory Standards, National

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AntiStreptolysin O

Anti-Streptolysin O (WHO)
The ASO assay's calibrators have been standardized against the WHO
(World Health Organization) NIBSC Anti-Streptolysin O standard material.

TYPE: Reagent
METHOD: Immunoturbidimetric
FORMAT: Liquid
LYNERITY: 20 - 1,000 IU/ML
SAMPLE: Serum

KIT COMPOSITION
R1: 4 X 20 mL buffer
R2: 4 X 20 mL Streptolysin O latex suspension
Tests / Kit: 400 (Hitachi 717)

REAGENT COMPOSITION
REAGENT #1 - BUFFER
Phosphate buffer saline (pH 7.43)
Polyethylene glycol 40 g / l
REAGENT #2 - LATEX
Glycine buffer (pH 8.2)
Latex sensitized to ASL 0.17%

PRINCIPLE
Anti-streptolysin O form immune complexes with sensitized latex.
Formed turbid suspension, the optical density is measured photometrically. The optical density depends on the content of the AS-O in the sample. The reaction produced by the "end-point" with a blank sample in the sample. The concentration is determined by the standard or non-linear calibration curve.

INTENDED USE:
For the quantitative determination of antibody to streptolysin O (ASO) in patient serum based on immunoturbidimetric assay as an aid in the diagnosis of Group A streptococcus infections.
FOR IN VITRO DIAGNOSTIC USE

REAGENT PREPARATION:
The reagents are ready to use and do not require reconstitution.

STORAGE AND HANDLING:
All reagents should be stored refrigerated (2-8°C). Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for one year from the date of manufacture, as indicated by the expiration date on the package and bottle labels. Opened reagents can be used for one month if stored at 2-8°C.

PERFORMANCE:
Precision
When a sample is assayed 20 times, the CV is ≤ 5%.

ACCURACY/CORRELATION:
A comparison of the AntiStreptolysin test and an AntiStreptolysin of another manufacturer gave such results:

y = 1.03x - 0.87
r = 0.994
n = 84
x = company A's ASO assay (IU/mL)
y = ASO (WHO) assay (IU/mL)

ASSAY RANGE:
20 - 1,000 IU/mL

INTERFERENCE:
Bilirubin F and C: No interference up to 40 mg/dL
Formazin: No interference up to 3,000 FTU
Hemoglobin: No interference up to 500 mg/dL

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KIT COMPOSITION:
R1: 2 X 22 mL buffer
R2: 1 X 7 mL antibody
Calibrator: 1 x 1 mL
Tests / Kit: 137 (Hitachi 717)
(Can Vary Based On Analyzer Model)

INTENDED USE:
For the quantitative determination of human Anti-Thrombin III (AT-III) in human plasma by immunoturbidimetric assay.
FOR IN VITRO DIAGNOSTIC USE.

PRINCIPLE:
The method of determining the activity of ATIII is based on the ability of ATIII to neutralize thrombin in the presence of heparin. ATIII activity in plasma is determined by adding a surplus of thrombin. This is accompanied by inhibition of thrombin ATIII-heparin complex is proportional to the amount of ATIII in the plasma. The remaining amount of thrombin catalyzes cleavage of para-nitroaniline (pNA) from synthetic chromogenic substrate. Absorption of free pNA, determined at 405 nm is inversely proportional to the activity ATIII.
The process goes as follows:
ATIII + heparin (excess) → ATIII-heparin
ATIII-heparin + thrombin (excess) → ATIII-heparin-thrombin + thrombin (residue)
Substrate pNA + thrombin (residue) → peptide + pNA

REAGENT PREPARATION:
The reagents are ready to use. They do not need to be reconstituted

REAGENT COMPOSITION:
R1
Lyophilized bovine thrombin with heparin. (6 IU / ml)
R2
SAR-PRO-ARG p-nitroanilide (3.1 mM)
R3
TAVS buffer (0.38 mol / l), pH 8.2;
EDTA (37.5 mmol / l);
Heparin (5 U / ml)
Stabilizer, sodium azide 0.09%

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PERFORMANCE:
Precision
The precision for the AT-III assay was determined using packaged reagents, pooled human plasma, and a Hitachi 717 analyzer.

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KIT COMPOSITION

R1 3 X 20 mL buffer
R2 1 X 20 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human apolipoprotein A-I, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of Apo A-I in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE

For the quantitative determination of human Apolipoprotein AI (Apo AI) in serum by immunoturbidimetric assay.
FOR IN VITRO DIAGNOSTIC USE.

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PERFORMANCE

Precision

The precision for the Apo AI assay was determined using packaged reagents, pooled human serum, and a Hitachi 704 analyzer.

Precision - Within Run

Sample I	Sample II	Sample III
N = 20	N = 20	N = 20
Mean = 83.1	Mean = 142.6	Mean = 189.0
S.D. = 1.6	S.D. = 2.3	S.D. = 2.9
C.V. = 1.87%	C.V. = 1.60%	C.V. = 1.51%

Precision - Between Runs

Sample I	Sample II	Sample III
N = 10	N = 10	N = 10
Mean = 85.9	Mean = 136.5	Mean = 185.9
S.D. = 1.1	S.D. = 1.437	S.D. = 2.9
C.V. = 1.28%	C.V. = 1.05%	C.V. = 1.57%

Accuracy / Correlation

A comparison of the Apo AI assay and a Sigma Apo AI Test Kit was performed using a Hitachi 704. The test results provided the following data:

$y = 0.960x + 4.776$
 $r = 0.970$
 $n = 55$
 $x = \text{Sigma Apo AI Test Kit}$
 $y = \text{Apo AI assay}$
 $x \text{ min} = 91$
 $x \text{ max} = 200$
 $x \text{ mean} = 146$
 $y \text{ min} = 89$
 $y \text{ max} = 213$
 $y \text{ mean} = 147$

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KIT COMPOSITION:
R1: 4 X 18 mL buffer
R2: 1 X 10 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human apolipoprotein A-II, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of Apo A-II in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE:

For the quantitative determination of human Apolipoprotein AII (Apo AII) in serum by immunoturbidimetric assay.
FOR IN VITRO DIAGNOSTIC USE.

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PERFORMANCE:

Precision:

When a sample that has an apolipoprotein AII concentration around 30 mg/dL is assayed 20 times (within-run), the absorbance C.V. is below 5%.

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KIT COMPOSITION 1:
R1: 3 X 20 mL buffer
R2: 1 X 20 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

KIT COMPOSITION 2:
Note: Large size kit, available by special order.
R1: 3 X 60 mL buffer
R2: 1 X 60 mL antibody
Tests / Kit: 600 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human apolipoprotein B, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of Apo B in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE:

For the quantitative determination of human Apolipoprotein B (Apo B) in serum by immunoturbidimetric assay.
FOR IN VITRO DIAGNOSTIC USE

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PERFORMANCE:

Precision

The precision for the Apo B assay was determined using packaged reagents, pooled human serum, and a Hitachi 717 analyzer.

Precision - Within Run

Sample I	Sample II	Sample III
N = 20	N = 20	N = 20
Mean = 36.0	Mean = 70.2	Mean = 120.2
S.D. = 0.7	S.D. = 1.5	S.D. = 1.7
C.V. = 2.06%	C.V. = 2.12%	C.V. = 1.44%

Precision - Between Runs

Sample I	Sample II	Sample III
N = 10	N = 10	N = 10
Mean = 58.6	Mean = 83.3	Mean = 124.5
S.D. = 1.3	S.D. = 1.9	S.D. = 1.4
C.V. = 2.25%	C.V. = 2.33%	C.V. = 1.09%

Accuracy / Correlation

A comparison of the Apo B assay and a Sigma Apo B Test Kit was performed using a Hitachi 704. The test results provided the following data:

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Apolipoprotein CII Reagent Set Аполипопротеин CII

Apolipoprotein CII

Type: Reagent
Method: Immunoturbidimetric
Format: Liquid
Linearity: 1 - 15 mg/dL
Sample: Serum

KIT COMPOSITION
R1: 3 X 20 mL buffer
R2: 1 X 20 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human apolipoprotein C-II, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of Apo C-II in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE:

For the quantitative determination of human Apolipoprotein CII (Apo CII) in serum by immunoturbidimetric assay

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PERFORMANCE:

Precision:

When a sample that has an apolipoprotein CII concentration around 4 mg/dL is assayed 20 times (within-run), the absorbance C.V. is below 5%.

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Precautions

For in vitro diagnostic use only.
Handle serum and blood specimens as potentially infectious samples and follow the guidelines established by the Centers for Disease Control (CDC), Atlanta, Georgia, for blood collection and handling.

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Apolipoprotein CIII Reagent Set Аполипопротеин CIII

Apolipoprotein C-III

Type: Reagent
Method: Immunoturbidimetric
Format: Liquid
Linearity: 3 - 30 mg/dL
Sample: Serum

KIT COMPOSITION:
R1: 3 X 20 mL buffer
R2: 1 X 20 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human apolipoprotein C-III, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of Apo C-III in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE:

For the quantitative determination of human Apolipoprotein CIII (Apo CIII) in serum by immunoturbidimetric assay.

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PERFORMANCE

Precision:

When a sample that has an apolipoprotein CIII concentration around 10 mg/dL is assayed 20 times (within-run), the absorbance C.V. is below 5%.

Assay Range:

3 - 30 mg/dL

INTERFERENCE:

Bilirubin: No interference up to 20 mg/dL
Hemoglobin: No interference up to 500 mg/dL
Lipemia: No interference up to 20%

Precautions

For *in vitro* diagnostic use only.
Handle serum and blood specimens as potentially infectious samples and follow the guidelines established by the Centers for Disease Control (CDC), Atlanta, Georgia, for blood collection and handling.

Materials Provided
Reagents.

Materials Required but not Provided
Biochemistry analyzer.

Analysis

See User Manual for your Biochemistry analyzer

Results registration

See User Manual for your Biochemistry analyzer

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KIT COMPOSITION:

R1: 3 X 20 mL buffer
R2: 1 X 20 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human apolipoprotein E, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of Apo E in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE:

For the quantitative determination of human Apolipoprotein E (Apo E) in serum by immunoturbidimetric assay.

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PERFORMANCE:

Precision

When a sample that has an apolipoprotein E concentration around 4 mg/dL is assayed 20 times (within-run), the absorbance C.V. is below 5%.

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INTENDED USE

For the quantitative determination of aspartate aminotransferase (AST) in human serum. The reagents are used in routine examination and monitoring of therapy and relapses.


INTRODUCTION

Serum aspartate aminotransferase (AST) also known as serum glutamic oxaloacetic transaminase (SGOT) is a tissue enzyme that catalyzes the exchange of amino and keto groups between alpha-amino acids and alpha-keto acids. AST is widely distributed in tissue principally cardiac, hepatic, muscle and kidney. Injury to these tissues results in the release of the AST (SGOT) enzyme to general circulation. Following a myocardial infarction, serum levels of AST (SGOT) are elevated and reach a peak 48 to 60 hours after onset. Hepatobiliary diseases, such as cirrhosis, metastatic carcinoma, and viral hepatitis also will increase serum AST levels.¹

The first kinetic assay of AST for diagnostic purposes was described by Karmen et al. in 1955, using a coupled reaction of malate dehydrogenase (MDH) and NADH.² This assay system was critically evaluated and optimized in 1960 by Henry et al.³ In 1977 the International Federation of Clinical Chemistry recommended a reference procedure for the measurement of AST activity based upon Karmen's procedures. The AST reagent applies the formulation recommended by the IFCC.

PRINCIPLE

The enzymatic reaction sequence employed in the assay of aspartate Aminotransferase is as follows: AST

AST catalyzes the transfer of an amino group between L-aspartate and 2-oxoglutarate. The oxalacetate formed in the first reaction is then reacted with NADH in the presence of malate dehydrogenase (MDH) to form NAD. AST activity is determined by measuring the rate of oxidation of NADH at 340 nm. Lactate dehydrogenase is included in the reagent to convert endogenous pyruvate in the sample to lactate during the lag phase prior to measurement.

REAGENT COMPOSITION

AST Liquid Reagents 1 and 2 come in separate containers, and both reagents are clear, colorless liquid in ready to use format. After combining AST Liquid R1 (Buffer Reagent) and AST Liquid R2 (Co-Enzyme) the reagent contains:

L-Aspartate 240 mmol/L
MDH (porcine muscle) >600 U/L
LDH (rabbit muscle) >600 U/L
Tris Buffer, pH 7.5 80 mmol/L
2-Oxoglutarate 12 mmol/L
NADH 0.18 mmol/L
Stabilizers and Preservatives

WARNINGS AND PRECAUTIONS

The reagents are for *in vitro* diagnostic use. Normal precautions exercised in handling laboratory reagents should be followed. The reagents contain sodium azide which may be toxic if ingested. Sodium azide may also react with lead and copper plumbing to form highly explosive metal azides. Refer to Material Safety Data Sheet for any updated risk, hazard or safety information.

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Example: 25 ml R1 + 5 ml R2

REAGENT STORAGE AND STABILITY

Reagents are stable until the expiration date on their respective labels, when properly stored at 2 - 8°C and protected from light. Reagents should appear clear and colorless. Discard if either appears cloudy or contains particulate matter. The working reagent is stable for 4 weeks at 2 - 8°C. The working

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reagent should be discarded if the initial absorbance, read against distilled water at 340 nm, is below 1.100.

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1. Spectrophotometer capable of absorbance reading at 340 nm and 1cm light path
2. Constant temperature block or bath, 37°C, or temperature controlled cuvette well
3. Accurate pipetting devices
4. Test tubes
5. Interval timer

SPECIMEN COLLECTION AND STORAGE

Non-hemolyzed serum is the specimen of choice, yet EDTA treated plasma or heparinized plasma can be used. Whenever possible, specimens should be separated and analyzed on the day of collection. Store serum in stoppered tubes. The enzyme in serum is reportedly stable for a minimum of 7 days at 2-8°C.

INTERFERING SUBSTANCES

Hemolysis must be avoided as the concentration of AST in red cells is roughly 10 times that of serum. Bilirubin levels up to 40 mg/dL and triglyceride levels up to 2000 mg/dL show no interference in this test. Certain drugs and other substances are also known to affect AST values.

MANUAL PROCEDURE

1. Prepare AST working reagent according to instructions.
2. Zero spectrophotometer at 340 nm with distilled water
3. For each sample and control, add 1.0 mL working reagent to cuvette or test tube and warm to 37°C for 3 minutes.
4. Add 100 µL (0.10 mL) serum to its respective tube and mix gently. Read and record absorbance at 1 minute. Continue incubating at 37°C and record absorbance again at 2 and 3 minutes. Rate should be constant.
5. Determine the average absorbance per minute (ΔA/min), multiply by factor 1768 for results in U/L.

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AUTOMATED PROCEDURE

Refer to appropriate application manual available

QUALITY CONTROL

It is recommended that control be included in each set of assays. Commercially available control material with established AST values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedural errors.

CALIBRATION

AST activity is based on the "micromolar extinction coefficient" of NADH at 340 nm (see "Results" section). The instrument manufacturer's calibration guidelines should be followed to calibrate your analyzer. Assaying the AST contents of a control serum with known AST values can be used to assure instrument calibration has been performed correctly.

RESULTS

Values are derived based on the "absorptivity micromolar extinction coefficient" of NADH at 340 nm (0.00622). Units per liter (U/L) of AST/GOT activity is that amount of enzyme which oxidizes one µmol/L of NADH per minute.

$$U/L = \frac{\Delta A / \text{Min}}{\text{Absorptivity}} \times \frac{\text{Total Volume}}{\text{Sample Volume}}$$

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LIMITATIONS

If the $\Delta A/\text{min}$ is greater than 0.342, dilute 1 part sample with 9 parts isotonic saline and re-assay. Multiply the result by 10. AST values for neonatal patients have not been established with this procedure. Grossly icteric or turbid specimen may require the use of a sample blank.

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PERFORMANCE CHARACTERISTICS

Comparison: A group of 125 sera ranging in AST activity from 13 - 399 U/L was assayed by the described AST method and by a similar commercially available AST reagent. Comparison of the results yielded a correlation coefficient of 0.999 and the regression equation was $y = 0.964x + 0.964$. (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T.)

Precision: Within-run precision was established by 30 assays on two different levels of commercial serum controls. Total Precision values were obtained by assaying the two commercial controls for 5 consecutive days.

	Within-Run	
	Serum 1	Serum 2
Mean AST (U/L)	41.7	115.2
Std. Deviation (U/L)	0.9	3.3
C.V. (%)	2.2	2.9

	Run to Run	
	Serum 1	Serum 2
Mean AST (U/L)	40.4	116.2
Std. Deviation (U/L)	0.87	4.88
C.V. (%)	2.14	4.19

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Linearity: Linear to 500 U/L at 37°C. Performed according to NCCLS Guideline EP6-P.

Sensitivity: Based on an instrument resolution of $A = 0.001$, the method presented shows a sensitivity of 2.65 U/L.

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REFERENCES

- Henry, J.B.: *Clinical Diagnosis and Management by Laboratory Methods*, W.B. Saunders and Co., Philadelphia, PA, p 361 (1974).
- Karmen A: A note on the spectrophotometric assay of glutamic-oxalacetic transaminase in human blood. *J Clin Invest* 34:131, 1955.
- Henry, R.J., et al.: *Am J Clin Path*, 34:381 (1960). The Committee on Enzymes of the Scandinavian Society for
- Clinical Chemistry and Clinical Physiology. *Scand J Clin Lab Invest* 32:291 (1974).
- Demetriou, JA et al. In *Clinical Chemistry - Principles and Techniques* 2nd ed. RJ Henry et al. Eds. Harper & Row, Hagerstown MD, 1974, p 873.
- International Federation of Clinical Chemistry. Provisional Recommendations on IFCC Methods for the Measurement of Catalytic Concentrations of

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Enzymes. Clin Chem 23: 887, 1977.

7. Young D.S. Effects of drugs on clinical laboratory tests. AACCC Press.
Washington D.C., 1990

8. Henry JB. Clinical Diagnosis and Management by Laboratory

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htimed.com

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High Technology, Inc

ACETAMINOPHEN

For the In Vitro Quantitative Determination of Acetaminophen in Serum

Summary and Principle

Acetaminophen (Tylenol, paracetamol, p-hydroxyacetanilide) is used in many formulations¹ as an analgesic generally with no adverse effects. Although the consequences of long-term treatment with acetaminophen are unclear, cases of hepatotoxicity and nephrotoxicity following chronic excessive use have been reported. However, in cases of overdose, it can cause severe hepatic damage leading to hepatic failure if untreated. Incidents of self poisoning as well as poisoning in small children have been increasing, now exceeding the incidence of aspirin poisoning. The management of acetaminophen overdose requires an early recognition of the drug in the blood stream at certain concentrations. The toxic range is generally reported as over 200 mg/L. N-acetylcysteine has been used as an antidote in conjunction with intensive supportive care. Early diagnosis of acetaminophen induced hepatotoxicity is important because initiation of therapy within 16 hours of ingestion lessens the potential for hepatic injury and decreases the mortality rate.² A specific, simple and quick method to determine cases of overdose is most helpful so that proper treatment can quickly be made available. An enzymatic method offers these advantages over other methods. Enzymes have been known to hydrolyze n-arylacylamides. In the Stanbio method, acetaminophen is hydrolyzed by an arylacylamidase enzyme to yield p-aminophenol and acetate. Subsequently, the p-aminophenol is measured colorimetrically by its conversion to an indophenol in the presence of o-cresol where periodate is used as a catalyst.

Reagents

Reagent A - Enzyme, Ref. No. 2401
Arylacylamidase < 7.5 KU/L
o-Cresol < 9 mmol/L
Reagent B - Catalyst, Ref. No. 2402
Sodium m-Periodate < 3.75 mmol/L
Standard - 300 mg/L, Ref. No. 2403
Acetaminophen 300 mg/L

Precautions

For in vitro diagnostic use only. Avoid skin contact with the reagents. If skin contact occurs, wash immediately with large amounts of water. Handle serum and blood specimens as potentially infectious samples and follow the guidelines established by the Centers for Disease Control (CDC), Atlanta, Georgia, for blood collection and handling.

Reagent Preparation:

Prepare the working reagent by adding 1 part reagent B to 2 parts reagent A, in that order. Mix before use. The working reagent is stable at 2-8 °C for 2 days or at 20-25°C for 4 hours when kept in stoppered amber vial.

Reagent Storage and Stability:

Reagents A and B are stable at 2-8°C for 24 months after manufacture when kept separated. The aqueous Acetaminophen Standard is stable at 2-8 °C for 24 months after manufacture. Do Not Freeze.

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Serum is the recommended specimen.
Sample Stability:

Interfering Substances: No interference was observed when the following drugs at a 1mM concentration were added to serum containing 100 mg/L (0.661 mM) acetaminophen.

acetylsalicylic acid
p-ethoxyacetanilide
amylorbarbital
pentazocine
caffeine
phenobarbital
chloridiazepoxide
promethazine
diazepam
salicylamide
diphenhydramine
salicylic acid
indomethacin

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p-aminosalicylic acid

Ethanol was tested at a concentration of 0.4% and found not to interfere. Acetylsalicylic acid and salicylic acid at concentrations of 5 mM were found not to interfere. N-acetyl-L-cysteine was tested at a concentration of 1000 mg/L and found not to significantly interfere. Concentrations of bilirubin at 10 mg/dL, creatinine at 25 mg/dL, or hemolytic serum with OD540 of 2.5 do not interfere significantly with the measurement of acetaminophen. However, at extremely high concentrations of bilirubin (25 mg/dL), the acetaminophen concentrations of 90, 180, and 270 mg/L, the observed values were 131 mg/L, 208 mg/L, and 284 mg/L.

EDTA at a concentration of 0.15% results in 10% lower reported values at concentrations of 100 and 200 mg/L acetaminophen. In serum contain 90 mg/L acetaminophen, a 1 mM concentration of the drugs amiflopridine, amphetamine, imipramine and p-ethoxyaniline gave biased values of -22%, -24%, -14% and +34%, respectively.

Procedure

The Acetaminophen LiquiColor® provided can be used in a 10 minute procedure as well as in a kinetic procedure, as described in Procedure 1 and Procedure 2, respectively. Both procedures give the same results provided that the spectrophotometer used reads absorbance to the third decimal place.

Procedure #1- 10 Minute

Test Conditions: Wavelength..... 615 nm

Temperature..... 37°C

Mode..... 10 minute rate

Sample Volume..... 100 µL

Reagent Volume..... 2.0 mL

Sample to reagent ratio..... 1/20

Procedure: 1) Prepare the required amount of working reagent by adding 1 part Reagent B to 2 parts Reagent A, in that order.

2) Incubate the working reagent at 37°C for 3 minutes

3) To two 3 mL cuvettes, add 2 mL of working reagent and measure OD615(T₀).

4) To cuvette #1, add 100 µL of test serum and measure OD615(T₁).

5) To cuvette #2, add 100 µL of acetaminophen standard and measure OD615(T₂).

6) Exactly 10 minutes after adding serum or standard, measure OD615(T₃).

7) Subtract T₀ from T₁ to get ΔOD10

Calculations:

ΔOD10 of test serum x 300 mg/L x dil of serum = Acetaminophen mg/L

ΔOD10 of standard

An acetaminophen concentration of 100 mg/L should give ΔOD10 of about

0.24 - 0.26 in the 10 minute test mode

PROCEDURE (for TC MATRIX)

Procedure #2 - Kinetic Test

Conditions: Wavelength..... 615nm

Temperature..... 37°C

Mode..... Kinetic

Sample Volume..... 100 µL

Reagent volume..... 2.0 mL

Sample to reagent ratio..... 1/20

Procedure: 1) Prepare the required amount of working reagent by adding 1 part Reagent B to 2 parts Reagent A, in that order.

2) Incubate the working reagent at 37°C for 3 minutes

3) To two 3 mL cuvettes, add 2 mL of working reagent.

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High Technology, Inc

- 5) To cuvette #2, add 100 µL of acetaminophen standard.
6) Determine the ΔOD/min between 3 and 5 mins at 37°C.

Calculation:

ΔOD/min of test serum x 300 mg/L x dil of serum = Acetaminophen mg/L
ΔOD/min of standard

NOTES: a) Precise measurements of temperature and time are required to obtain accurate results.

b) The working reagent when first prepared has a light yellow color which gradually deepens over the lifetime of the working reagent. Test results are not affected by this coloration.

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Expected Values

Although toxic manifestations have been observed at serum concentrations of > 100 mg/L. The therapeutic range varies, and has been reported at 10-30 mg/L.

Performance Characteristics

Linearity: Both procedures described above are linear to 300 mg/L acetaminophen. For higher concentrations, dilute the sample with distilled water. Repeat the assay and multiply the results by the dilution factor.

Sensitivity: Concentrations of acetaminophen of 90, 100 and 110 mg/L can be clearly distinguished at the 99.9% confidence limit with either method. Due to the variation between laboratory instruments, and the reduced absorbance change with acetaminophen concentrations below 10 mg/L, test results below 10 mg/L should be reported as "<10 mg/L."

Precision: Precision studies were conducted using two serum pools containing 50 and 200 mg/L acetaminophen. The following results are an average of eighteen determinations.

10 Minute Test: Standard
Mean (mg/L) Deviation CV(%)
Within Day 50.0 0.71 1.42
202.9 3.05 1.50
Between Day 50.0 1.79 3.58
202.7 6.36 3.14

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INTENDED USE

For the quantitative determination of total protein concentration in human serum.

INTRODUCTION

Serum protein is involved in the maintenance of normal distribution of water between blood and tissues through osmotic pressure. Low protein is primarily caused by malnutrition, impaired synthesis, loss (as by hemorrhage), or excessive protein catabolism. Elevated protein levels are caused mainly by dehydration.¹

The determination of total protein in serum makes use of the Biuret color reaction, known since 1878. Past attempts to stabilize the cupric ions in the alkaline reagent were unsuccessful until the addition of sodium potassium tartrate as a complexing agent.² The present method for quantitative determination of total protein in serum is based on the method proposed by the American Association for Clinical Chemistry³ (AACC) and National Committee for Clinical Laboratory Standards (NCCLS).⁴

PRINCIPLE

The enzymatic reaction sequence employed in the assay of Total Protein is as follows:

Protein in serum forms a blue colored complex when reacted with cupric ions in an alkaline solution. The intensity of the violet color is proportional to the amount of protein present when compared to a solution with known protein concentration.

REAGENT COMPOSITION

1. Total Protein Reagent: Sodium Hydroxide 600 mM, Copper Sulfate 12 mM, Sodium Potassium Tartrate 32 mM, Potassium Iodide 30 mM, and non-reactive ingredients.

2. Total Protein Standard: Bovine Albumin Ft. V with preservative 5.0 g/dl.

REAGENT STORAGE AND STABILITY

Store Total protein reagent at room temperature (18

Total protein standard at refrigerator (2 - 8°C).

REAGENT DETERIORATION

The reagent should be discarded if there is turbidity, or the presence of a black precipitate, which indicates reagent deterioration. The reagent should be a clear, pale blue solution.

WARNINGS AND PRECAUTIONS

1.

For in vitro diagnostic use

CAUTION: In vitro diagnostic reagents may be hazardous. Handle in accordance with good laboratory procedures which dictate avoiding ingestion, and eye or skin contact.

Specimens should be considered infectious and handled appropriately.

Avoid ingestion. DO NOT PIPETTE BY MOUTH.

The reagent contains sodium hydroxide that is corrosive. In case of contact with skin, flush with water. For eyes, seek medical attention.

attention.

SPECIMEN COLLECTION

1. Test specimens should be serum and free from hemolysis.
2. Gross hemolysis will cause elevated results because of the released hemoglobin as well as the increase in background color.
3. Lipemic sera cause elevated results and should be run with a serum blank.
 - a. Place 1.0 ml 0.9% saline in test tube.
 - b. Add 0.02 ml (20 ul) sample
 - c. Zero spectrophotometer with 0.9% saline.
 - d. Read and record absorbance of serum blank.
 - e. Subtract blank absorbance from test absorbance.
 - f. Calculate as usual.
4. Samples with bromsulphthalein (BSP) will result in falsely elevated results.⁵

Protein in serum is stable for one (1) week at room temperature (15 - 30°C) and for at least one (1) month refrigerated (2 - 8°C) when guarded against evaporation.

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1. Accurate pipetting devices (3 ml, 50 ul)
2. Timer
3. Test tube and rack
4. Spectrophotometer

GENERAL INSTRUCTIONS

The reagent for Total Protein is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

AUTOMATED PROCEDURE

Refer to appropriate application manual available.

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1. Label tubes as Blank, Standard, Control, Patient, etc.
2. Pipette 3.0 ml of reagent into each tube.
3. Add 0.05 ml (50 ul) of standard and patients to appropriate tubes and mix by inversion.
4. Let the tubes stand at room temperature (15 - 30°C) for ten (10) minutes.
5. Set spectrophotometer at 540 nm and zero instrument with the reagent blank. (Wavelength range 500-560 nm)
6. Read and record absorbance of each tube.

* TC MULTI-PURPOSE CALIBRATOR MAYBE USED TO REPLACE STANDARD.

NOTE: Final color is stable for sixty-minutes at room temperature.
ALTERNATE VOLUMES: 20 ul (0.02 ml) sample to 1.2 ml reagent. Calculations remain the same.

LIMITATIONS

1. The reagent is linear to 15.0 g/dl. Samples with values above 15.0 g/dl should be diluted 1:1 with 0.9% saline, re-run, and the result multiplied by two (2).

2. The Buret procedure is not sensitive at low ranges (< 1 g/dl).

3. Not for use with urine or spinal fluid specimens.

CALCULATIONS

Abs. Of Unknown x Conc. Of Standard = Total Protein (g/dl) Abs. Of Standard

Example:

Absorbance of unknown = 0.350

Absorbance of standard = 0.400

Concentration of standard = 5 g/dl

$0.350 \times 5 = 4.38$ g/dl

0.400

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established total protein values may be routinely used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction, or procedural errors.

EXPECTED VALUES

6.2-8.5 g/dl

1. The effect of posture, when blood is drawn, varies with the individual but recurrent values are usually lower than ambulatory. Differences may be as much as 1.2 g/dl

2. It is strongly recommended that each laboratory establish its own range of expected values.

PERFORMANCE CHARACTERISTICS

1. Linearity: 1.0-15.0 g/dl
2. Comparison: A comparison study when performed between this procedure and another procedure based on the same principle resulted in a correlation coefficient of 0.95 with a regression equation of $y = 0.86 + 1.02$.

Mean (g/dl)	Within Run	
	S.D.	C.V. %
6.8	0.12	1.8
3.7	0.08	2.1

Mean (g/dl)	Run-to-Run	
	S.D.	C.V. %
6.8	0.17	2.4
3.7	0.14	3.7

REFERENCES

1. Peters, T. and Biamonte, G.T., *Selected Methods for the Small Clinical Chemistry Laboratory*. Faubler, W.R., and Melles, S., Ed.
2. Gomali, A. et al. *J. Clin. Chem* 177:75(1949).
3. Doumas, B.T., et al. *Clin. Chem* 27:1642 (1981)
4. NCCLS Approved Standards: ACS-1, *Specification for Standardized Protein Solution (Bovine Serum Albumin)*. 2nd ed., National Committee for Clinical Laboratory Standards. 771 E Lancaster Ave., Villanova, PA 19 - 85, (1979).
5. Henry, R.J., et al. *Clinical Chemistry Principles and Techniques*. Harper & Row, N.Y., 415 (1974).
6. Young, D.S., et al. *Clin. Chem* 21:10-432, (1975).
7. Tietz, N.W., *Fundamentals of Clinical Chemistry*. W.B. Saunders, Philadelphia, PA 298, (1976).

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Intended use

Reagent for measuring β 2-microglobulin by latex-enhanced turbidimetric immunoassay. Application to both Serum/Plasma and Urine.

Method

Beta2-microglobulin, contained in patient samples, reacts to the coated latex particles, causing their agglutination. The degree of the agglutination can be determined by optical density at 546 nm. Increased turbidity of the solution is directly proportional to the concentration of beta 2-microglobulin in the sample.

Reagents

R1. Buffer for analysis

Tris buffer, 20 mmol/L, pH 8.2

BSA, 5 g/l

Sodium azide, 0.95 g/l

R2. Reagent with antibodies

Latex particles coated with β 2-M antibodies

Dynamic measuring range: 0.2 – 80 mg/L (serum), 30–8000 μ g/L (urine).

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INTENDED USE

For the quantitative determination of total bilirubin in human serum.

INTRODUCTION

Bilirubin is a metabolite of the heme portion of heme proteins, mainly hemoglobin. Bilirubin is excreted into the intestine and bile from the liver. The site of the catabolism of hemoglobin is the reticuloendothelial system (RES). Bilirubin is then released into the bloodstream where it binds tightly to albumin and is transported to the liver. Upon uptake by the liver, bilirubin is conjugated with glucuronic acid to form bilirubin mono and diglucuronide, which are water-soluble metabolites. The metabolites are then excreted in the bile.¹

Elevation of total serum bilirubin may occur due to (1) excessive hemolysis or destruction of the red blood cells (e.g. hemolytic disease of the newborn), (2) liver diseases (e.g. hepatitis and cirrhosis), (3) obstruction of the biliary tract (e.g. gallstones).¹

Most chemical methods for the determination of total bilirubin are based on the reaction between diazotized sulfanilic acid and bilirubin to produce azobilirubin, which absorbs maximally at 560 nm. Such tests are often run in the presence or absence of an organic solvent e.g., methanol to distinguish free bilirubin from conjugated bilirubin on a differential solubility basis.

PRINCIPLE

Bilirubin reacts with diazotized sulfanilic acid to produce azobilirubin, which has an absorbance maximum at 560 nm in the dimethyl sulfoxide (DMSO) solvent. The intensity of the color produced is directly proportional to the amount of total bilirubin concentration present in the sample.

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REAGENT PREPARATION

Total bilirubin working reagent: Add 0.05 ml (50 µl) of sodium nitrite reagent per 1.0 ml of total bilirubin reagent and mix reagent.

Example: 0.5ml sodium nitrite/0.05ml total bilirubin reagent, 1 ml sodium nitrite/20ml total bilirubin reagent, etc.

REAGENT STORAGE

1. Store reagent at refrigerator temperature (2-8°C).
2. Combined working reagent can be stored for up to eight (8) hours when kept in an amber bottle at room temperature.
3. Do not freeze reagents.
4. Avoid exposure to direct sunlight.

REAGENT DETERIORATION

- The reagent should be discarded if:
1. Sodium Nitrite reagent has a yellow discoloration.
 2. Working reagent fails to achieve assigned assay values of fresh control sera.

SPECIMEN COLLECTION AND STORAGE

1. Hemolysis interferes with the test, i.e. Hemolyzed samples should be avoided since they may give falsely low values.³

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- All specimens for this assay must be carefully protected from light.¹
- Bilirubin in serum is stable for 4-7 days when stored in the dark at 2-8 °C.²

INTERFERENCES

- Young et al. gives an exhaustive list of drugs and other substances known to affect the circulating level of bilirubin.³
- In this assay, as in all laboratory procedures, materials, which come in contact with specimens, should be clean and free of contamination by heavy metals, detergents, and other chemicals.

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- Test tubes, rack
- Pipettes
- Timers
- Appropriate automated chemistry analyzer or spectrophotometer capable of measuring at 560 nm

TOTAL BILIRUBIN PROCEDURE (AUTOMATED)

Refer to appropriate instrument application instructions.

TOTAL BILIRUBIN PROCEDURE (MANUAL)

- Label test tubes, "Blank, Standard, Control, Patient". Each tube requires a blank tube.
- Dispense 1.0ml of total bilirubin reagent to all blank tubes.
- Prepare a working reagent. See "Reagent Preparation."
- Dispense 1.0ml of the working reagent into the labeled test tubes, except blank tubes.
- Add 0.1ml (100 µl) of each standard, control and sample to its respective tube. Mix well.
- Allow all tubes to stand for five (5) minutes at room temperature.
- Set the wavelength of the instrument at 560nm. Zero with reagent blank. (Wavelength range: 500-550nm).
- Read and record absorbance of all tubes.

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PROCEDURE NOTES

- The final color produced is stable for 60 minutes.
- For pediatric samples with bilirubins over 3.0 mg/dl, use 0.05ml (50 µl) of sample and then multiply the result by two (2).
- If the spectrophotometer being used requires a final volume greater than 1.0ml for accurate readings, 3ml of reagent and 0.2ml (200 µl) of sample could be used.

CALCULATIONS

Abs. = absorbance

Total bilirubin (mg/dl):

$$\frac{\text{Abs. of unknown} - \text{Abs. of blank}}{\text{Abs. of calibrator} - \text{Abs. of calibrator blank}} \times \text{Conc. of calibrator}$$

Example:

Unknown	= 0.565
Absorbance of unknown blank	= 0.024
Absorbance of calibrator	= 0.480
Absorbance of calibrator blank	= 0.004
Concentration of calibrator	= 5.0 mg/dl

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INTRODUCTION

Bilirubin is a metabolite of the heme portion of heme proteins, mainly hemoglobin. Bilirubin is excreted into the intestine and bile from the liver. The site of the catabolism of hemoglobin is the reticuloendothelial system (RES). Bilirubin is then released into the bloodstream where it binds tightly to albumin and is transported to the liver. Upon uptake by the liver, bilirubin is conjugated with glucuronic acid to form bilirubin mono and diglucuronide, which are water-soluble metabolites. The metabolites are then excreted in the bile.¹ Elevation of total serum bilirubin may occur due to (1) excessive hemolysis or destruction of the red blood cells (e.g. hemolytic disease of the newborn), (2) liver diseases (e.g. hepatitis and cirrhosis), (3) obstruction of the biliary tract (e.g. gallstones).¹

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REAGENT DETERIORATION

The reagent should be discarded if:

1. Sodium Nitrite reagent has a yellow discoloration.
2. Working reagent fails to achieve assigned assay values of fresh control sera.

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SPECIMEN COLLECTION AND STORAGE

1. Hemolysis interferes with the test, i.e. Hemolyzed samples should be avoided since they may give falsely low values.
2. All specimens for this assay must be carefully protected from light.
3. Bilirubin in serum is stable for 4-7 days when stored in the dark at 2-8 °C.

INTERFERENCES

1. Young et al. gives an exhaustive list of drugs and other substances known to affect the circulating level of bilirubin.
2. In this assay, as in all laboratory procedures, materials, which come in contact with specimens, should be clean and free of contamination by heavy metals, detergents, and other chemicals.

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1. Test tubes, rack
2. Pipettes
3. Timers
4. Appropriate automated chemistry analyzer or spectrophotometer capable of measuring at 560 nm.

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TOTAL BILIRUBIN PROCEDURE (MANUAL)

1. Label test tubes, "Blank, Standard, Control, Patient". Each tube requires a blank tube.
2. Dispense 1.0ml of total bilirubin reagent to all blank tubes.
3. Prepare a working reagent. See "Reagent Preparation".
4. Dispense 1.0ml of the working reagent into the labeled test tubes, except blank tubes.
5. Add 0.1 ml (100 µl) of each standard, control and sample to its respective tube. Mix well.
6. Allow all tubes to stand for five (5) minutes at room temperature.
7. Set the wavelength of the instrument at 560nm. Zero with reagent blank. (Wavelength range: 500-550nm).
8. Read and record absorbance of all tubes.

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PROCEDURE NOTES

1. The final color produced is stable for 60 minutes.
2. For pediatric samples with bilirubins over 3.0 mg/dl, use 0.05ml (50 µl) of sample and then multiply the result by two (2).
3. If the spectrophotometer being used requires a final volume greater than 1.0ml for accurate readings, 3ml of reagent and 0.2ml (200 µl) of sample could be used.

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$\frac{\text{Abs. of unknown} - \text{Abs. of blank}}{\text{Abs. of calibrator} - \text{Abs. of calibrator blank}} \times \text{Conc. of calibrator}$

Example:

Unknown = 0.565
Absorbance of unknown blank = 0.024
Absorbance of calibrator = 0.480
Absorbance of calibrator blank = 0.004
Concentration of calibrator = 5.0 mg/dl

Then:
 $\frac{0.565 - 0.024}{0.480 - 0.004} \times 5 = \frac{0.541}{0.476} \times 5 = 5.68 \text{ mg/dl}$

PROCEDURE LIMITATIONS

1. Serum with values above 20mg/dl must be diluted to 1:1 with isotonic saline, reassayed and the final result multiplied by two (2).

2. Serum hemoglobin levels of up to 1.0g/dl do not interfere with results.

QUALITY CONTROL

It is recommended that commercially available control material with known total bilirubin values be included in each set of assays.

EXPECTED VALUES

Infants (after one month) and adults 0.2 - 1.0 mg/dl.

PERFORMANCE

1. Linearity: 20 mg/dl

2. Sensitivity: Based on an instrument resolution of A = 0.001, this procedure has a sensitivity of 0.01 mg/dl

3. Comparison: Studies between the present method and similar method yielded a correlation of R² = 0.99 and a regression equation of y = 1.13x - 0.09.

4. Precision:

Run-to-run: Two commercial control sera were assayed for a period of 30 days and the following day-to-day precision was obtained.

	Level I	Level II
Mean (mg/dl), N=22	0.68	6.41
S.D.	0.05	0.63
C.V.	7.3%	9.7%

Within Run: Two commercial control sera were assayed 20 times and the following within run precision was obtained.

	Level I	Level II
Mean (mg/dl)	0.67	6.28
S.D.	0.02	0.15
C.V.	2.9%	2.3%

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High Technology, Inc

INTENDED USE

For the quantitative determination of direct (conjugated) bilirubin concentration in serum or plasma on Auto Analyzer analyzers.

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An increase in the formation or retention of Bilirubin in the body results

in increased levels of serum Bilirubin and jaundice. This

hyperbilirubinemia is classified as either pre-hepatic, hepatic or

post-hepatic depending on the principal cause of condition. Therefore,

determination of the total Bilirubin and its conjugated (direct) Bilirubin

is important for the differential diagnosis of hyperbilirubinemia and

Bilirubin reacts with diazotized sulfanilic acid to produce azobilirubin,

which has an absorbance maximum at 560 nm in the aqueous solution.

The intensity of the color produced is directly proportional to the

amount of direct bilirubin concentration present in the sample.

The Auto Analyzer System automatically proportions the appropriate

sample and reagent volumes into the cuvette. The system monitors the

change in absorbance at 560 nanometers. This change in absorbance is

directly proportional to the concentration of direct bilirubin in the

sample and is used by the Auto Analyzer System to calculate and express

direct bilirubin concentration.

Direct Bilirubin + diazo + H⁺ → azobilirubin (red color)

REAGENT CONTENTS:

Each kit contains: Four Direct Bilirubin Reagent 1 (4*30 ml)

Two Direct Bilirubin Reagent 2 (2*6 ml)

Instruction Insert.

REAGENT PREPARATION

No preparation is required.

REAGENT COMPOSITION

Sulfanilic acid: 27 mM;

Sodium nitrite: 0.12 mM;

HCl: 51 mM.

Also non-reactive chemicals for optimal system performance.

REAGENT STORAGE AND STABILITY

Direct Bilirubin Reagent stored unopened at 2°C to 8°C is stable until

the expiration date shown on the bottle label. Once opened Direct Bilirubin

Reagent is stable for 30 days, unless the expiration date

is exceeded. Bilirubin Calibrator stored at 2°C to 8°C is stable until the

expiration date shown on the bottle label.

DO NOT FREEZE.

SPECIMEN COLLECTION AND HANDLING

1. The test can be performed on serum, plasma. For serum, blood is drawn into a tube which does not contain anticoagulant and allow clotting. The serum is then separated from the clot. A maximum limit of two hours from the time of collection is recommended.

2. Separated serum or plasma should not remain at room temperature longer than 8 hours. If assays are not completed within 8 hours, serum and plasma should be stored at 2°C to 8°C. If assays are not completed within 48 hours, or the separated sample is to be stored beyond 48 hours, samples should be frozen at -15°C to -20°C. Frozen samples should be thawed only once. Analyte deterioration may occur in samples that are repeatedly frozen and thawed.

3. For plasma, add whole blood directly into a tube containing anticoagulant.

Acceptable anticoagulants are listed in the "LIMITATIONS" section.

CALIBRATION

1. Calibrator required: MULTI Calibrator.

2. The system must have a valid calibration in memory before controls or patient samples can be run.

3. The Auto Analyzer system will automatically perform checks on the calibration and produce data at the end of calibration.

Note: Refer to the Auto Analyzer manual for further instructions on calibrating the instrument.

MATERIALS NEEDED BUT NOT SUPPLIED WITH

REAGENT KIT

MULTI Calibrator

At least two levels of control material.

LIMITATIONS

1. The anticoagulants Potassium Oxalate, Sodium Fluoride EDTA, Lithium Heparin, Sodium Heparin, Ammonium Heparin and Sodium Citrate were found to be compatible with this method.

INTERFERENCE

1. Hemoglobin may interfere with this methodology.

2. Lipemic samples >3+ should be ultra-centrifuged and the analysis performed on the infranate.

3. On this method, refer to the work of Young for a review of drug and comprehensive list of substances effect on level.

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High Technology, Inc

PRECAUTIONS:

1. For in vitro diagnostic use only.
2. Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
3. Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests

PROCEDURES

TEST NAME:	DBIL	R1:	310
TEST NO.		R2:	30
FULL NAME:	Direct Bilirubin	SAMPLE VOLUME:	10
REFERENCE NO.		R1 BLANK:	/
ANALY TYPE:	Endpoint	MIX REAG. BLANK:	0.1 - 1.3
PR. WAVE:	546 nm	CONCENTRATION:	/
SECON. WAVE:	/	LINEARITY LIMIT:	0 - 20
TREND:	Ascending	SUBSTRATE LIMIT:	/
REACT. TIME:	0 - 15	FACTOR:	/
INCUBATE TIME:	3	PROZONE CHECK:	/
UNIT:	mg/dl	Q1 / Q2 / Q3 / Q4 /	
PRECISION:	0.1	PC / ABS /	

Calibration Type: Calibrate+ Reag.Blank Calibration Rule: Two-point linear

Analytical Range: 0.1-10 mg/dl

For Direct Bilirubin analyte by Direct Bilirubin Reagent on Four TC Matrix System, this method has been demonstrated to be linear from 0.1-10 mg/dl.

Accuracy: Comparison study was performed on Auto Analyzer System from 40 samples. Beckman Coulter Direct Bilirubin Reagent was used to compare with Direct Bilirubin Reagent. The results of this study in yield a correlation coefficient of 0.99 with a regression equation of $y=0.99x-0.01$

Precision: Within Run precision for Direct Bilirubin Reagent Set was determined following a modification of NCCLS EP5-A Two commercial human serum, Direct Bilirubin Calibrator were assayed on Auto Analyzer System for 25 times.

Sample	Sample 1	Sample 2
N	25	25
Mean (mg/dl)	0.3	1.4
Standard Deviation (mg/dl)	0.01	0.02
Coefficient of Variation (%)	4.5	4.0

Run-Day precision for Direct Bilirubin Reagent was determined following a modification of NCCLS EP5-A Two commercial human serum, Direct Bilirubin Calibrator were assayed on Auto Analyzer Systems five times per day for five days for the total of 25 values.

Sample	Sample 1	Sample 2
N	25	25
Mean (mg/dl)	0.3	1.4
Standard Deviation (mg/dl)	0.01	0.02
Coefficient of Variation (%)	4.5	4.0

REFERENCES:

1. Malloy, H.T. and Evelyn, K.A. 1960, J.Biol. Chem. 119, 48.
2. Valdez, O.S. et al. 1971 J. Pediatr. 79, 1015.
3. Tietz, N.W., "Specimen Collection and Processing: Sources of Biological Variation," Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders, Philadelphia, PA (1994)
4. National Committee for Clinical Laboratory Standards. Approved Guideline, NCCLS publication C28-A, Villanova, PA (1994).
5. Henry, J. B., ed., Clinical Diagnostics and Management by Laboratory Methods, 18th Edition, W.B. Saunders, Philadelphia.
6. Tietz, N.W., ed., Clinical Guide to Laboratory Tests, 2nd Edition, W.B. Saunders, Philadelphia, PA (1990)
7. National Committee for Clinical Laboratory Standards, Method Comparison and Bias

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1. Malloy, H.T. and Evelyn, K.A. 1960, J.Biol. Chem. 119, 48.
2. Valdez, O.S. et al. 1971 J. Pediatr. 79, 1015.
3. Tietz, N.W., "Specimen Collection and Processing: Sources of Biological Variation," Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders, Philadelphia, PA (1994)
4. National Committee for Clinical Laboratory Standards. Approved Guideline, NCCLS publication C28-A, Villanova, PA (1994).
5. Henry, J. B., ed., Clinical Diagnostics and Management by Laboratory Methods, 18th Edition, W.B. Saunders, Philadelphia.
6. Tietz, N.W., ed., Clinical Guide to Laboratory Tests, 2nd Edition, W.B. Saunders, Philadelphia, PA (1990)
7. National Committee for Clinical Laboratory Standards, Method Comparison and Bias Estimation Using Patient Samples. Tentative Guideline, NCCLS Publication EP9-T,

High Technology, Inc

Estimation Using Patient Samples; Tentative Guideline, NCCLS Publication EP9-T, Villanova, PA (1993)
8. National Committee for Clinical Laboratory Standards, Precision Performance of Clinical Chemistry Devices; Tentative Guideline, 2nd Edition, NCCLS publication EP5-T2, Villanova, PA (1992)

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Villanova, PA (1993)
8. National Committee for Clinical Laboratory Standards, Precision Performance of Clinical Chemistry Devices; Tentative Guideline, 2nd Edition, NCCLS publication EP5-T2, Villanova, PA (1992)

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INTENDED USE

For the quantitative determination of direct bilirubin in human serum. For in vitro diagnostic use.

INTRODUCTION

Bilirubin is a metabolite of the heme portion of heme proteins, mainly hemoglobin. Normally it is excreted into the intestine and bile from the liver. The site of the catabolism of hemoglobin is the reticuloendothelial system (RES). Bilirubin is then released into the bloodstream where it binds tightly to albumin and is transported to the liver. Upon uptake by the liver, bilirubin is conjugated with glucuronic acid to form bilirubin mono and diglucuronide that are water soluble metabolites. The metabolites will react with aqueous diazo reagent and are commonly referred to as "direct bilirubin".¹

Elevation of total serum bilirubin may occur due to (1) excessive hemolysis or destruction of the red blood cells e.g. hemolytic disease of the newborn, (2) liver diseases e.g. hepatitis and cirrhosis (3) obstruction of the biliary tract e.g., gallstones.¹ There is information in the literature indicating elevated levels of direct bilirubin in patients with liver or biliary tract diseases: even though, total bilirubin levels are normal.² Therefore, the greatest diagnostic value of direct bilirubin assays stem from their ability to indicate occult liver disease. Most chemical methods for the determination of total bilirubin are based on the reaction between diazotized sulfanilic acid and bilirubin to produce azobilirubin, which absorbs maximally at 560 nm. Such tests are often run in the presence and absence of an organic solvent e.g., methanol to distinguish free bilirubin from conjugated bilirubin on a differential solubility basis.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. Specimens should be considered infectious and handled appropriately.
3. Do not pipette reagents by mouth. Avoid contact reagent with eyes, skin and clothing. Do not ingest. Wash hands after use.

PRINCIPLE

Bilirubin reacts with diazotized sulfanilic acid to produce azobilirubin, which has an absorbance maximum at 560 nm in the aqueous solution. The intensity of the color produced is directly proportional to the amount of direct bilirubin concentration present in the sample.

REAGENTS

1. Bilirubin Reagent: Sulfanilic Acid 32mM, Hydrochloric Acid 165mM.
2. Bilirubin Nitrite Reagent: Sodium Nitrite 60mM.
3. Bilirubin Calibrator: N-1-naphthyl ethylenediamine dihydrochloride salt. (5 mg/dl).

REAGENT PREPARATION

Direct bilirubin working reagent: Add 0.005 ml (5 µl) of sodium nitrite reagent per 1.0 ml of direct bilirubin reagent and mix reagent. Example: 0.05ml sodium nitrite/0.05ml direct bilirubin reagent, 0.1ml sodium nitrite/20ml direct bilirubin reagent, etc.

REAGENT STORAGE

1. All reagents are stored at room temperature (15 - 30°C).
2. Combined working reagent can be stored for up to eight (8) hours when kept in an amber bottle at room temperature.
3. Do not freeze reagents.
4. Avoid exposure to direct sunlight.

REAGENT DETERIORATION

1. The reagent should be discarded if:
2. Sodium Nitrite reagent has a yellow discoloration.
3. Working reagent fails to achieve assigned assay values of fresh control sera.

SPECIMEN COLLECTION AND STORAGE

1. Hemolysis interferes with the test, i.e. Hemolyzed samples should be avoided since they may give falsely low values.¹
2. All specimens for this assay must be carefully protected from light.¹
3. Bilirubin in serum is stable for 4-7 days when stored in the dark at 2-8 °C.

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INTERFERENCES

1. Young et al. give an exhaustive list of drugs and other substances known to affect the circulating level of bilirubin.
2. In this assay, as in all laboratory procedures, materials, which come in contact with specimens, should be clean and free of contamination by heavy metals, detergents, and other chemicals.

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1. Cuvettes
2. Pipettes
3. Timers
4. Appropriate automated chemistry analyzer or spectrophotometer capable of measuring at 560 nm.

DIRECT BILIRUBIN PROCEDURE (AUTOMATED)

Refer to appropriate instrument application instructions

DIRECT BILIRUBIN PROCEDURE (MANUAL)

1. Label test tubes, "Blank, Standard, Control, Patient" Each tube requires a blank tube.
 2. Dispense 1.0ml of direct bilirubin reagent to all blank tubes.
 3. Prepare a working reagent. See "REAGENT PREPARATION"
 4. Dispense 1.0ml of the working reagent into the labeled test tubes.
 5. Add 0.1ml (100µl) of each standard, control, and sample to its respective tube. Mix well.
 6. Allow all tubes to stand for five (5) minutes at room temperature.
 7. Set the wavelength of the instrument at 560nm. Zero with reagent blank. (Wavelength range: 500-560).
 8. Read and record absorbance of all tubes.
- * HTI - MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

PROCEDURE NOTE

- The final color produced is stable for 60 minutes.
- For pediatric samples with bilirubin over 3.0 mg/dl, use 0.05ml (50 µl) of sample and then multiply the result by two (2).
2. If the spectrophotometer being used requires a final volume greater than 1.0ml for accurate readings, 3ml of reagent and 0.2ml (200 µl) of sample may be used.

CALCULATIONS

Abs = absorbance

$$\frac{\text{Abs. of unknown} - \text{Abs. of unknown blank}}{\text{Abs. of calibrator} - \text{Abs. of calibrator blank}} \times \text{Concentration of calibrator (mg/dl)} = \text{Direct Bilirubin (mg/dl)}$$

Example:	
Absorbance of unknown	= 0.132
Absorbance of unknown blank	= 0.120
Absorbance of calibrator	= 0.450
Absorbance of calibrator blank	= 0.000
Concentration of calibrator	= 5.0 mg/dl

$$\frac{0.132 - 0.120}{0.450 - 0.000} \times 5 = \frac{0.012}{0.450} \times 5 = 0.13 \text{ mg/dl}$$

PROCEDURE LIMITATIONS

1. Sera with values above 20 mg/dl must be diluted 1:1 with isotonic saline, reassayed and the final answer multiplied by two (2).
2. Serum hemoglobin levels of up to 1.0 g/dl do not interfere with results.

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QUALITY CONTROL

Normal and abnormal control sera of known concentrations of direct bilirubin should be analyzed routinely with each group of unknown specimens.

EXPECTED VALUES^{1,6}

Infant (after one month) and adults: 0.0-0.5 mg/dl. It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE

- Linearity:** 20 mg/dl
- Sensitivity:** Based on an instrument resolution of 0.001 absorbance, the present procedure has a sensitivity of 0.01 mg/dl.
- Comparison:** A comparison study between the present method with an available commercial product using the same identical method on forty (40) fresh serum samples and two commercial serum controls, ranging from 0.10 mg/dl to 0.67 mg/dl yielded a coefficient of 0.98 and a regression equation of $Y = 0.99X + 0.01$.

4. Precision studies:

Day-to-Day Precision: Two commercial control sera were assayed for a period of 21 days and the following day-to-day precision was obtained.

	Level 1	Level 2
Mean (mg/dl) N = 16	0.13	0.33
S.D.	0.01	0.04
C.V.%	7%	12%

Within Run Precision: Two commercial control sera were assayed 20 times and the following within run precision was obtained.

	Level 1	Level 2
Mean (mg/dl) N = 20	0.10	0.23
S.D.	0.01	0.03
C.V.	10.0%	13.0%

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Intended Use


For the quantitative determination of β -hydroxybutyrate in serum or plasma. For *in vitro* diagnostic use only.

Summary and Explanation of Test

Ketosis is a common feature in acutely ill patients. In subjects suffering from starvation, acute alcohol abuse, or diabetes mellitus ketosis can result in severe life threatening metabolic acidosis. The presence and degree of ketosis can be determined by measuring blood levels of β -hydroxybutyrate. Ordinarily, β -hydroxybutyrate is the ketoacid present in the greatest amount in serum. It accounts for approximately 75% of the ketone bodies that also contain acetoacetate and acetone. During periods of ketosis, β -hydroxybutyrate increases even more than the other two ketoacids, acetoacetate and acetone, and has been shown to be a better index of ketoacidosis including the detection of subclinical ketosis. In diabetics, the measurement of β -hydroxybutyrate as well as blood glucose is needed for the assessment of the severity of diabetic coma and is essential for the exclusion of hyperosmolar non-ketotic diabetic coma. Moreover, the insulin requirements are often based on the extent of the existing hyperketonemia shown by the blood levels of β -hydroxybutyrate. A specific enzymatic assay for β -hydroxybutyrate is therefore extremely important in the assessment of ketosis.

Principle

Enzymatic quantitation of β -hydroxybutyrate by β -hydroxybutyrate dehydrogenase has been reported. In the High Technology Inc. method, β -hydroxybutyrate (D-3-hydroxybutyrate) in the presence of NAD gets converted to acetoacetate and NADH at a pH 8.5 by β -hydroxybutyrate dehydrogenase (D-3-hydroxybutyrate dehydrogenase). At this pH, the reaction is favored to the right. The NADH produced is converted to color using INT and diaphorase.

Reagents

The β -hydroxybutyrate reagent kit contains the following.

1. 1 x 50ml R1 containing β -hydroxybutyrate dehydrogenase and diaphorase enzymes.
2. 1 x 8ml R2 containing NAD, INT and oxalate.
3. 1 x 3ml standard containing 1mM sodium D-3-hydroxybutyrate

Reagent Storage and Stability

All of the above reagents are stable stored at 2-8°C until the expiration date stated on the labels. Do not use the reagents past their expiration date. Alterations in the physical appearance of the reagents or values of control materials outside of the manufacturer's acceptable range may be an indication of reagent instability.

Reagent Preparation

R1 and R2 reagents are supplied as ready to use liquids. To prepare a single liquid working reagent mix 10 parts R1 with 1.5 parts R2 (ex. 10.0ml R1 with 1.5 ml R2).

Warnings and Precautions

For *in vitro* diagnostic use only. Avoid skin contact with reagents. If this occurs, wash immediately with water.

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Serum or plasma collected with EDTA, heparin or sodium fluoride can be used in the assay. Serum or plasma β -hydroxybutyrate levels are stable at least one week if kept refrigerated. (2-8°C)

Interferences

Lactate dehydrogenase and lactate have been shown to interfere with the assay. The incorporation of oxalic acid in this reagent eliminates this interference as reported. No significant changes in values were observed when the following analytes were added to serum containing 0.5mM β -hydroxybutyrate.

Analyte % Recovery
Glucose (2000 mg/dl) 96
Acetoacetic acid (5 mM) 96
Creatinine (5 mg/dl) 106

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In addition, hemolyzed serum with an OD540nm of 2.0 was added to the test and found not to interfere.

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1. 10 ml Micropipettor
2. 1 ml Pipette or Dispenser
3. Glass or plastic Test Tubes to hold 1.0 ml
4. β -hydroxybutyrate controls
5. Temperature controlled incubator
6. Spectrophotometer capable of reading at 505nm

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Procedure (Manual)

The High Technology Inc. enzymatic β -hydroxybutyrate reagent provided can be used in a five-minute procedure utilizing a spectrophotometer that reads at 505 nm.

Conditions: Wavelength 505nm

Temperature 37°C

Mode Endpoint

Reagent Volume 1.0 ml

Sample volume 25 ul

Incubation 5 min. Prepare working reagent according to "Reagent Preparation" instructions.

2. Label tubes, "Standard", "Control", "Patient", etc.

3. Pipette 1.0 ml of working reagent into required number of labeled test tubes. Incubate at 37°C for 5 minutes.

4. Add 25 ul Standard, controls and patients to the appropriate tests and incubate at 37°C for 5 minutes.

5. Zero spectrophotometer with a reagent blank.

6. Read and record absorbances of all tubes at 505 nm.

NOTES:

a) Precise measurement of temperature, wavelength, and time are required to obtain accurate results.

b) The test can also be performed at 25°C for 10 minutes using the same procedure as above.

Calibration

Use the β -hydroxybutyrate standard provided.

Calculations

β -hydroxybutyrate (mM) = OD (5min) serum x Standard concentration OD (5min) standard

To convert results to mg/dL, divide the value obtained in mM by 0.096

Limitations

1. See "Interferences"

2. The procedure described above is linear to 4.5mM (46.8 mg/dl) β -hydroxybutyrate. For higher concentrations dilute the sample with deionized water, repeat the assay and multiply the results by the dilution factor.

Quality Control

The reliability of test results should be monitored whenever patient samples are assayed using a standard and quality control materials analyzed in the same manner employed for the unknowns. We suggest the use of commercially available β -hydroxybutyrate controls with an assayed range. If controls do not fall into the assayed range, patient values from that run should not be reported. The run should be repeated, making sure that all mixing and handling instructions are strictly followed.

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Expected Values

The quantitation of β -hydroxybutyrate is important in cases of ketoacidosis.

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In studies of healthy individuals who had fasted for 12 hours before blood collection, the range of β -hydroxybutyrate was found to be from 0.02 mM (0.2mg/dl) to 0.27mM (2.81mg/dl). 4,5 Other ranges have also been reported.

Performance

1. Linearity: The β -hydroxybutyrate assay is linear to 4.5 mM (46.8 mg/dl).
2. Comparison: A study of 57 samples using High Technology Inc. and the Sigma method gave the following correlation: (y =High Technology Inc., x =Sigma); $y = 1.05x + 0.25$ mM, $N = 57$, $r = 0.9922$, $S_{xy} = 0.1055$ mM.
3. Precision: Studies were conducted using two serum pools containing 0.25mM (2.6 mg/dl) and 1.0mM (10.4 mg/dl) β -hydroxybutyrate. The following results are averages of eighteen determinations:

Within day:

Level	Mean (mM)	Std. Dev. (mM)	% C.V.
1	0.29	0.005	1.7
2	1.09	0.015	1.4

Day to day:

Level	Mean	Std. Dev.	% C.V.
1	0.26	0.014	5.2
2	1.05	0.018	1.7

4. Sensitivity: Concentrations of β -hydroxybutyrate of 0.18, 0.28, and 0.38 mM (1.8, 2.9, and 3.9 mg/dl) can be clearly distinguished at the 99% confidence limit.

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KIT COMPOSITION

R1: 3 X 18 mL buffer
R2: 2 X 7 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human haptoglobin, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of haptoglobin in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE

For the quantitative determination of human haptoglobin in human serum by immunoturbidimetric assay. Haptoglobin is a protein that binds hemoglobin, the oxygen-carrying pigment in red blood cells. Measurement of haptoglobin may aid in the diagnosis of hemolytic diseases (diseases in which the red blood cells rupture and release hemoglobin) relating to the formation of hemoglobin-haptoglobin complexes and certain kidney diseases.
FOR IN VITRO DIAGNOSTIC USE.

REAGENT PREPARATION

The reagents are ready to use. They do not need to be reconstituted.

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PERFORMANCE

Precision

The precision was determined using packaged reagents, pooled human serum, and a Hitachi 717 analyzer.

Precision - Within Run

Sample I	Sample II
N = 20	N = 20
Mean = 68.15	Mean = 366
S.D. = 0.88	S.D. = 5.5
C.V. = 1.3%	C.V. = 1.5%

Precision - Between Runs

Sample I	Sample II
N = 8	N = 8
Mean = 85.9	Mean = 225.6
S.D. = 0.8	S.D. = 1.4
C.V. = 0.972%	C.V. = 0.624%

Accuracy / Correlation

A comparison of Haptoglobin assay and an Incstar Haptoglobin Test Kit was performed using a Hitachi 717. The test results provided the following data:

$y = 0.937x - 9.25$			
$r = 0.993$			
$n = 60$			
$x = \text{Incstar Haptoglobin Test Kit}$			
$y = \text{K-ASSAY}^\circledR \text{ Haptoglobin assay}$			
$x \text{ min} =$	34	$y \text{ min} =$	4
$\text{max} =$	266	$\text{max} =$	267
$\text{mean} =$	109	$\text{mean} =$	93

Assay Range

20 - 220 mg/dL

INTERFERENCE

Bilirubin F and C :	No interference up to 20 mg/dL
Hemoglobin :	No interference up to 500 mg/dL
Lipemia :	No interference up to 20%

Precautions

For in vitro diagnostic use only.
Handle serum and blood specimens as potentially

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INTENDED USE

For the Kinetic Quantitative Determination of γ -Glutamyl Transferase in Serum for Manual and/or Automated Procedures

SUMMARY AND PRINCIPLE

Gamma-glutamyl transferase (γ -GT) is one of a large group of enzymes known as peptidases. Although renal tissue has the highest level of γ -GT, the major source of the enzyme present in serum is of hepatic origin. Elevated levels of γ -GT are found in association with hepatobiliary and pancreatic disorders, alcoholism and heavy drinkers, in myocardial disorders and in diabetics.

Unlike alkaline phosphatase activity, the serum γ -GT activity remains normal in diseases affecting bone and during normal bone growth. Therefore, a rise in serum γ -GT activity may be considered as a sensitive and more specific indicator of liver disease than alkaline phosphatase activity.

The HTI γ -GT procedure has been optimized according to Szasz.

γ -GT

L- γ -glutamyl-3-carboxy-4-nitroanilide + glycylglycine $\rightarrow$
L- γ -glutamylglycylglycine + 5-amino-2-nitrobenzoate
 γ -GT catalyzes the transfer of a γ -glutamyl group from L- γ -glutamyl-3-carboxy-4-nitroanilide. The rate of liberation of 5-amino-2-nitrobenzoate is directly related to the γ -GT activity in the sample and is quantitated by measuring the increase in absorbance at 405 nm.

REAGENTS

γ -GT Buffer (R1)

Tris, pH 8.25 100 mmol/L
Glycylglycine 100 mmol/L

γ -GT Substrate (R2)

L- γ -Glutamyl-3-Carboxy-4-Nitroanilide 4.0 mmol/L

PRECAUTIONS

The reagents are for "In Vitro Diagnostic Use". Normal precautions exercised in handling laboratory reagents should be followed. The reagents contain sodium azide that may be toxic if ingested. Sodium azide may also react with lead and copper plumbing to form highly explosive metal azides. Refer to Material Safety Data Sheet for any updated risk, hazard or safety information.

REAGENT PREPARATION

Buffer and Substrate liquid reagents are supplied ready-to-use. Prepare Working Reagent in the ratio of 5 parts Buffer (R1) to 1 part Substrate (R2), (i.e., 25 mL Buffer and 5 mL Substrate).

REAGENT STORAGE AND STABILITY

Reagents are stable until the expiration date on their respective labels, when properly stored at 2-8°C and protected from light. R1 should appear clear/colorless while R2 should appear clear/yellow. Discard if either appears cloudy or contains particulate matter. The Working Reagent is stable for 4 weeks at 2-8°C or 5 days at room temperature (15-25°C). The Working Reagent should be discarded if the initial absorbance, read against distilled water at 405 nm, is above 0.800.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Spectrophotometer capable of absorbance reading at 405 nm and 1 cm light path
2. Constant temperature block/bath, 37°C, or temperature controlled cuvette
3. Accurate pipetting devices
4. Test tubes
5. Interval timer

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Complexing anticoagulants such as citrate, oxalate, fluoride and must be avoided since they inhibit y-GT activity². The loss of y-GT activity is minimized by storing the samples refrigerated for up to 7 days or frozen up to 2 months³. Bilirubin levels up to 40 mg/dL and triglyceride levels up to 2000 mg/dL show no interference in this test.

INTERFERING SUBSTANCES

y-GT is an inducible enzyme. Consequently patients who are receiving antiepileptic drugs or aminopyrine show elevated y-GT activity. Chronic use of ethanol also increases serum y-GT activity⁴⁻⁵. Certain drugs and other substances are also known to affect y-GT values.

MANUAL PROCEDURE

1. Prepare y-GT Working Reagent according to instructions.
2. Zero spectrophotometer at 405 nm with distilled water.
3. For each sample and control, add 1.0 mL Working Reagent to cuvette or test tube and warm to 37°C for 3 minutes.
4. Add 100 µL (0.10 mL) serum to its respective tube and mix gently.
5. Read and record absorbance at 1 minute. Continue incubating at 37°C and record absorbance again at 2 and 3 minutes. Rate should be constant.
6. Determine the average absorbance per minute (ΔA/min), multiply by factor 1158 for results in U/L.

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AUTOMATED PROCEDURE

Special adaptations for automated analyzers are available by contacting HTI's Customer Service Department.

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established y-GT values may be used for quality control. The assigned value of the control material must be confirmed by this methodology. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction or procedural errors.

CALIBRATION

y-GT activity is based on the "micromolar extinction coefficient" of 5-amino-2-nitrobenzoate at 405 nm (see "Results" section). The instrument manufacturer's calibration guidelines should be followed to calibrate your analyzer.

RESULTS

Values are derived based on the "absorptivity micromolar extinction coefficient" of 5-amino-2-nitrobenzoate at 405 nm (0.0095). Units per liter (U/L) y-GT activity are the amount of enzyme that produces one nmol/L of 5-amino-2-nitrobenzoate per minute.

$$U/L = \frac{\Delta A/\text{min}}{\text{Absorptivity}} \times \frac{\text{Total Volume}}{\text{Sample Volume}}$$

$$U/L = \frac{\Delta A/\text{min}}{0.0095} \times \frac{1.10}{0.10}$$

$$U/L = \Delta A/\text{min} \times 1158$$

LIMITATIONS

If the ΔA/min. is greater than 0.259, dilute 1 part sample with 5 parts isotonic saline and re-assay. Multiply results by 6.

EXPECTED VALUES

Normal Range:

Males: 0 - 50 U/L (37°C)

Females: 0 - 30 U/L (37°C)

This range should serve only as a guideline. It is recommended that

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each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

PERFORMANCE CHARACTERISTICS

Comparison: A group of 63 sera ranging in γ -GT activity from 15 - 714 U/L was assayed by the described γ -GT method and by a similar commercially available γ -GT reagent. Comparison of the results yielded a correlation coefficient of 1.000 and the regression equation was $y = 0.94x + 1.9$. (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T.)

Precision: Within-run precision was established by 20 assays on three different levels of commercial serum controls. Total Precision values were obtained by assaying the 3 commercial controls for 5 consecutive days.

Within-Run

	Serum 1	Serum 2	Serum 3
Mean γ -GT (U/L)	40	74	206
Std. Deviation (U/L)	1.0	0.9	1.3
C.V. (%)	2.5	1.2	0.6

Total Precision

	Serum 1	Serum 2	Serum 3
Mean γ -GT (U/L)	42	72	204
Std. Deviation (U/L)	0.6	0.6	0.7
C.V. (%)	1.5	0.9	0.4

Precision studies were performed according to NCCLS Tentative Guideline, EP5-T.

Linearity: Linear to 300 U/L at 37°C. Performed according to NCCLS Guideline EP6-P.

Sensitivity: Based on an instrument resolution of $A = 0.001$, the method presented shows a sensitivity of 1.0 U/L.

REFERENCES

1. Tietz, N.W. Fundamentals of Clinical Chemistry, 2nd ed. W.B. Saunders Co., Philadelphia, PA, p. 622 (1976)
2. Whitfield, J.B., Moss DW, Neale, G, Orme, M and Breckenridge, A. Brit Med J 1, 316, 1973
3. Rosalki SB, Rau D, Lehman D, Prentice M: Determination of γ -glutamyltranspeptidase activity and its clinical applications. Ann Clin Biochem 7, 143, 1970
4. Rosalki SB, Tarlow D, Rau D: Plasma γ -glutamyltranspeptidase elevation in patients receiving enzyme inducing drugs. Lancet 2: 376, 1971
5. Bartels H, Hauck W, Vogel I: aminopyrine an effective modifier of liver and serum γ -glutamyltranspeptidase. J Pediatr 86: 298, 1975
6. Young DS, Effects of drugs on clinical laboratory tests. AACCC Press, Washington D.C. 1960
7. Tietz Textbook of Clinical Chemistry, 2nd ed. W.B. Saunders Co., Philadelphia, PA, p. 851 (1994)

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Hemoglobin Reagent Set Гемоглобин

INTENDED USE

Hemoglobin reagent set is used for the quantitative determination of hemoglobin in human blood.

INTRODUCTION

Hemoglobin is a porphyrin-iron (II) protein compound that transports oxygen from the lungs to body tissues where it is utilized for energy metabolism. Measurements of hemoglobin from venous or capillary blood aid in the detection of a variety of conditions which alter the normal hemoglobin concentration of blood, e.g. anemia or polycythemia. The determination of iron content in whole blood is the most accurate method for assessing blood hemoglobin. Of the various methods used, cyanmethemoglobin is the most widely accepted. It is this internationally adapted method that is employed in this procedure. 1

PRINCIPLE

In the cyanmethemoglobin method, erythrocytes are lysed by a stomatolytic agent in the presence of a surfactant and release their hemoglobin into solution. Hemoglobin is oxidized to methemoglobin by ferricyanide, and the methemoglobin is converted into the stable cyanmethemoglobin by addition of KCN. The absorbance of cyanmethemoglobin is measured at 540 nm and color intensity is proportional to hemoglobin concentration.

REAGENT COMPOSITION

When reconstituted as directed, the reagent for hemoglobin contains the following:

1. Hemoglobin reagent: Potassium ferricyanide 0.5 mM, potassium cyanide 0.7 mM, buffers and stabilizers included.
2. Standard: Methemoglobin (80 mg/dl) dissolved in cyanmethemoglobin reagent. This amount is equivalent to 15.0 g/dl hemoglobin. This standard has been referenced against a CAP (College of American Pathologists) certified standard to its concentration and further checked by using the known molar absorptivity of cyanmethemoglobin.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
- CAUTION:** *In vitro* diagnostic reagents may be hazardous. Handle in accordance with good laboratory procedures which dictate avoiding ingestion, and eye or skin contact.
2. Contains cyanide. Poison - may be fatal if swallowed. DO NOT PIPETTE BY MOUTH.
3. Do not mix with acids. Discarding with large volumes of water.
4. Specimens should be considered infectious and handled appropriately.
5. Use distilled or deionized water where indicated.

REAGENT PREPARATION

Reagent comes in a ready to use form.

REAGENT STORAGE AND STABILITY

Store the hemoglobin reagent and standard at room temperature (15 - 30°C).

REAGENT DETERIORATION

1. Do not use hemoglobin reagent if:
1. It has become a different color than yellow.
2. The reagent becomes turbid or a precipitation forms.

SPECIMEN COLLECTION

1. Use whole blood with EDTA as an anticoagulant.
 2. Oxalate, citrate or heparin may also be used as anticoagulants.
 3. Capillary or venous blood may be collected if used before clotting occurs.
 4. Whole blood mixed well with an anticoagulant appears stable for one (1) week at room temperature (15 - 30°C).
- INTERFERING SUBSTANCES**
1. Substances that cause turbidity will falsely elevate the hemoglobin value. These include lipids, abnormal plasma proteins (macroglobulinemia) or erythrocyte stroma.
 2. A review by Young *et al.* reveals the numerous drugs that exert an

in vivo effect to decrease blood hemoglobin.3

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1. Accurate pipetting devices
2. Timer
3. Test tubes/rack
4. Spectrophotometer with ability to read at 540 nm

GENERAL INSTRUCTIONS

The reagent for Hemoglobin is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

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MANUAL PROCEDURE

1. Dispense 2.0 ml of hemoglobin reagent into test tubes labeled "blank", "control", "patient", etc.
2. Place 0.01 ml (10 µl) of sample into respective tubes. Mix.
3. Allow all tubes to stand for three (3) minutes at room temperature.
4. To a tube labeled standard, place 2.0 ml of standard.
5. Set spectrophotometer to 540 nm and zero with the reagent blank. (Wavelength range: 520 - 550 nm).
6. Read and record absorbance values of all tubes.
7. See CALCULATIONS to obtain values.

NOTES:

1. For spectrophotometers requiring greater volumes for proper reading, use 4.0 ml reagent and 0.02 ml (20 µl) sample. Follow above instructions.
2. Final color appears quite stable but should be read within one (1) hour to avoid evaporation.

LIMITATIONS

1. This procedure measures hemoglobin and its derivatives except sulfhemoglobin.
2. Specimens with values above 20.0 g/dl must be re-run using one half the sample volume. Multiply final results by two (2).

CALIBRATION

Use hemoglobin standard provided.

CALCULATIONS

Abs. = Absorbance

$\frac{\text{Abs. of unknown}}{\text{Abs. of standard}} \times \text{Cone. of Standard (g/dl)} = \text{Value (g/dl)}$

Abs. of standard

Example: If a 15 g/dl standard has an absorbance of 0.602 and the absorbance of the unknown is 0.480 then:
 $\frac{0.480}{0.602} \times 15.0 = 11.9 \text{ g/dl}$
 0.602

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established hemoglobin values may be routinely used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedural errors.

EXPECTED VALUES

Adult Males	13.0-18.0 g/dl
Adult Females	11.0-16.0 g/dl
Children	10.0-14.0 g/dl
Newborns	14.0 - 23.0 g/dl

Factors such as age, race, exercise, season and altitude are reported to influence the values of normal ranges. The above range should serve

only as a guideline. Each laboratory should establish its own range

PERFORMANCE CHARACTERISTICS

1. Linearity: 20 g/dl

2. Sensitivity: Based on an instrument resolution of 0.001 absorbance, the present procedure has a sensitivity of 0.03 g/dl.

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3. Comparison: Studies conducted against a similar procedure yielded a coefficient of correlation of 0.98 with a regression equation of $y = 1.03x - 0.48$ on samples with values from 7.2 to 17.9 g/dl (n = 20).

4. Precision:

Within Run: Two samples of human blood were assayed twenty (20) times and the following within run precision was obtained.

	Mean (g/dl)	S.D.	C.V.%
Normal	13.8	0.6	4.6
Abnormal	10.2	0.3	3.4

Run-to-Run: Two samples of human blood were assayed for five (5) consecutive days and the following run to run precision was obtained.

	Mean (g/dl)	S.D.	C.V.%
Normal	14.3	12.3	0.5
Abnormal	12.3	0.5	4.3

REFERENCES

1. Eilers, R.J., *Am. J. Clin. Pathol.* 47:212 (1967)
2. Tietz, N.W., *Fundamentals of Clinical Chemistry*, 2nd ed. Saunders Co., Philadelphia p411 (1976)
3. Youm, D.S. et al., *Clin. Chem.* 21:10, (1975)
4. Henry, R.F. et al., *Principles and Techniques in Clinical Chemistry*, 2nd ed., Harper and Row, Hagerstown, MD p. 1135(1974)
5. Wolf, P.L., *Practical Clinical Hematology*, John Wiley and Sons, NY, p. 144 (1973)

Manufactured by:
HTI MEDICAL
109 PRODUCTION RD
WALPOLE, MA 02081
U.S.A.

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High Technology, Inc

Intended Use

For the quantitative determination of Glycohemoglobin (HbA1c) in blood by cation exchange resin. The test is to be used to monitor long-term glucose control in diabetes mellitus.

Summary and Explanation of Test

Throughout the circulatory life of the red cell, glycohemoglobin is formed continuously by the adduction of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al showed glycohemoglobin in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that glycohemoglobin serve as an indicator of metabolic control of the diabetic, since glycohemoglobin levels approach normal values for diabetics in metabolic control 2,3,4. Glycohemoglobin has been defined operationally as the "fast fraction" hemoglobins (HbA1a, A1b, A1c) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA0. The present glycohemoglobin procedure employs a weak binding cation-exchange resin for the rapid separation of glycohemoglobin (fast fraction) from non-glycosylated hemoglobin. Over 80% of the labile fraction of glycohemoglobin is removed during the separation step in this procedure due to the inclusion of the borate buffer system.⁵

Principle

A hemolyzed preparation of the whole blood is mixed continuously for five minutes with a weak binding cation-exchange resin. During this time, HbA0 binds to the resin. After the mixing period, a filter is used to separate the supernatant containing the glycohemoglobin from the resin. (Note: This binding is temperature dependent. Therefore, a standard should be included in each run.) The percent glycohemoglobin is determined by measuring the absorbance at 415 nm (405-420 nm acceptable) of the glycohemoglobin fraction and the total hemoglobin fraction. The ratio of the two absorbances gives the percent glycohemoglobin.

Reagents

40 Test Kit containing:

- 1 x 120ml Bottle 8mg/ml Cation-exchange Resin in a borate buffer, pH 6.9.
- 1 x 30ml Bottle Glycohemoglobin Lysing Resin, 10mM Potassium Cyanide, surfactant added.
- 40 serum separators.

Reagent Storage

Store reagent at room temperature (21-26°C).

Expiration Dating

All reagents are stable to expiration date stated on the labels. Do not use the reagents past their expiration date.

Reagent Deterioration

Alterations in the physical appearance of the reagents or values of control materials outside of the manufacturer's acceptable range may be an indication of reagent instability.

Instruments

Use a spectrophotometer able to read at 415nm with linearity to at least 1:5 O.D. units. (405-420 nm is acceptable.)

Precautions

- 1. This reagent is for *in vitro* diagnostic use only.
- 2. Not for internal or external use in humans or animals.
- 3. Lysing reagent contains cyanide (poison). Do not ingest. Do not mix with acid or HCN gas may be released.

Specimen Collection and Preparation

Special preparation of the patient is unnecessary. Fasting specimens are not required. No special additives or preservatives other than anticoagulants are required. Collect venous blood with EDTA using aseptic technique. All human specimens should be regarded as potentially biohazardous. Therefore, universal precautions should be used in specimen handling (gloves, lab garments, avoid aerosol production, etc.).

Storage

Glycohemoglobin in whole blood collected with EDTA is stable for one week at 2-8°C.

Interferences

Samples that are severely lipemic may cause elevated results. It has been reported that bilirubinemia may interfere with ion-exchange methods.⁶ Fetal Hemoglobin (HbF) may interfere in this assay. Blood samples with total hemoglobin greater than 18 g/dl should be diluted x 2 with deionized water before assay.

Materials Provided

Refer to "Reagents"

Materials Required but Not Provided

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High Technology, Inc

2. 500 µl, 3 ml and 5 ml Pipettes or Dispensers
3. 13x100 mm Glass Tubes
4. Glass or plastic Test Tubes to hold 0.6 ml and 5 ml
5. Rocker or Rotator (e.g. Miles Inc., Diagnostics Division 4651)
6. Glycohemoglobin standard (Cat. No. G7540-STD)
7. Glycohemoglobin controls (Cat. No. G7540-2)
8. Deionized water

Procedure

a. Hemolysate Preparation

1. Dispense 500 µl Lysing Reagent into tubes labeled: Standard, Control, etc. Note: Plastic or glass tubes of appropriate size are acceptable.
2. Place 100 µl of the well-mixed blood sample, standard or control into the appropriately labeled lysing reagent tube. Mix.
3. Allow to stand for 5 minutes or until complete lysis is evident.

b. Glycohemoglobin Preparation

1. Dispense 3.0 ml of Glycohemoglobin Cation-exchange Resin into 13 x 100 mm glass tubes labeled: Standard, Control, etc. Note: Before use, mix the resin by inverting at least 10 times. Swirl the bottle after addition to each 5 tubes.
2. Add 100 µl of the hemolysate (from step A3) to resin reagent.
3. Position the filter separators in the tubes so that the rubber sleeve is approximately 1 cm above the liquid level.
4. Place the tubes on the rocker or rotator and mix continuously for 5 minutes.
5. Remove the tubes from the rocker or rotator.
6. Push the filter separators into the tubes until the resin is firmly packed.
7. The supernatant may be poured into another tube or directly into a cuvette for absorbance measurement.
8. Zero the instrument at 415 nm (405-420 nm acceptable) with deionized water as the blank.
9. Read and record the absorbance values for Standard, Control, etc. These readings are for glycohemoglobin.
- c. Total Hemoglobin Fraction
1. Dispense 5.0 ml deionized water into plastic, or glass tubes labeled: Standard, Control, etc.
2. Place 20 µl of the hemolysate (from step A3) into the appropriately labeled tube of total hemoglobin diluent. Mix.
3. Adjust the instrument to zero absorbance at 415 nm (405-420 nm acceptable) with deionized water as the blank.
4. Read and record the absorbance values for Standard, Control etc. These readings are for total hemoglobin.

NOTE: This glycohemoglobin assay should be performed at room temperature, 21-26°C. The final reaction products for glycohemoglobin and total hemoglobin appear quite stable. However, the test samples should be read within an hour before evaporation becomes significant.

Limitations

1. This assay should not be used for the diagnosis of diabetes mellitus.
2. This method can be influenced by temperature. Patient specimens should always be assayed with a calibrator included in the run to eliminate temperature influences.
3. Glycosylated HbS and HbC bind more tightly than HbA, and produce lower values. Other hemoglobinopathies (e.g., beta-thalassemia and hemolytic anemia also produce lowered results.)
4. Results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, uremia (carbamylated Hb), alcoholism, ingest large doses of aspirin 7,8,9,10

Quality Control

The reliability of test results should be monitored whenever patient samples are assayed using a standard and quality control materials analyzed in the same manner employed for the unknowns. We suggest the use of commercially available glycohemoglobin controls with an assayed range. If controls do not fall into the assayed range patient values from that run should not be reported. The run should be repeated, making sure that all mixing and handling instructions are strictly followed.

Linearity of the assay should be verified with a commercial linearity check set, or dilutions of a high specimen, at least every six months.

Calculations

Results for the unknowns and controls are calculated as follows:

For each sample, calculate the ratio (R) of the glycohemoglobin absorbance to the total hemoglobin absorbance. Use the following equation to determine unknown concentrations:

$$\text{Unknown (\%)} = R (\text{Unk}) \times \text{Std Conc. (\%)} / R (\text{Std})$$

Example: A standard containing 10.0% glycohemoglobin had Abs. = 0.490 for the glycohemoglobin fraction and Abs. = 0.560 for the total hemoglobin fraction. An unknown sample had glycohemoglobin Abs. = 0.750 and total hemoglobin Abs. = 0.625. The glycohemoglobin concentration of the unknown is:

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Expected Values

Normal: 6.0 - 8.3%

The normal range represents the 95% confidence interval for 75 subjects with normal glucose values and no history of diabetes. Each laboratory should establish its own expected values. In using glycated Hb to monitor diabetic patients, results should be interpreted individually. That is, the patient should be monitored against him or herself. There is a 3-4 week time lag before % glycohemoglobin reflects changes in blood glucose level.

Performance

1. Linearity: The glycohemoglobin assay shows linearity for glycohemoglobin levels in the range of 4.0%-20%.

2. Comparison: A study on normal and abnormal human specimens between this glycohemoglobin procedure and a widely used column procedure yielded a correlation coefficient of 0.970 and a linear regression equation of $y=1.10x-3.00$ ($n=36$, range=5.9-14.2%).

3. Precision: Within Run: The intra assay precision was established by assaying blood with normal and elevated glycohemoglobin levels twenty times each. Level Mean Std. Dev. % C.V.

.V. Normal 7.7 0.24 3.1 Elevated

12.8 0.21 1.6

Run To Run: The inter run precision was established by assaying blood with normal and elevated glycohemoglobin levels for ten runs conducted over a five day period.

Level Mean Std. Dev. % C.V.

Normal 8.0 0.32 4.0

Elevated 14.8 0.45 3.0

4. Sensitivity: This glycohemoglobin procedure has a sensitivity of 0.02% glycohemoglobin per 0.001 units of absorbance.

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Intended Use

For the quantitative determination of Hemoglobin A1c (HbA1c) in human blood. The determination of HbA1c is most commonly performed for the evaluation of glycemic control in diabetes mellitus. HbA1c values provide an indication of glucose levels over the preceding 4-8 weeks. A higher HbA1c value indicates poorer glycemic control. For *in vitro* diagnostic use only.

Summary and Explanation of Test

Throughout the circulatory life of the red cell, Hemoglobin A1c is formed continuously by the addition of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivett et al. showed Hemoglobin A1c in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that Hemoglobin A1c serve as an indicator of metabolic control of the diabetic, since Hemoglobin A1c levels approach normal values for diabetics in metabolic control. HbA1c has been defined operationally as the last fraction* hemoglobins (HbA_{1a}, HbA_{1b}, HbA_{1c}) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA₀. The present procedure utilizes an antigen and antibody reaction to directly determine the concentration of the HbA1c.

Principle

This method utilizes the interaction of antigen and antibody to directly determine the HbA1c in whole blood. Total hemoglobin and HbA1c have the same unspecific absorption rate to latex particles. When mouse anti-human HbA1c monoclonal antibody is added (R2), latex-HbA1c-mouse anti human HbA1c antibody complex is formed. Agglutination is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA1c absorbed on to the surface of latex particles. The amount of agglutination is measured as absorbance. The HbA1c value is obtained from a calibration curve.

Reagents

R1: Latex 0.13%, Buffer, stabilizer.

R2 (When combined): Mouse anti-human HbA1c monoclonal antibody 0.05mg/ml, goat anti-mouse IgG polyclonal antibody 0.08mg/dl, Buffer, stabilizers.

Hemolysis reagent: water and stabilizers. (Included in 40ml kit, not 120ml kit)

Reagent Storage

Store all reagents refrigerated at 2-8° C.

Reagent Preparation

R1 and Hemolysis reagents are supplied as ready to use liquids. R2 is prepared by pouring the entire contents of the R2b vial into the R2a vial. Mix gently.

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Instruments

Refer to specific instrument application for suggested settings.

Precautions

1. This reagent is for *in vitro* diagnostic use only.

2. Not for internal or external use in humans or animals.

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Special preparation of the patient is unnecessary. Fasting specimens are not required. No special additives or preservatives other than anticoagulants are required. Collect venous blood with EDTA using aseptic technique. All human specimens should be regarded as potentially biohazardous. Therefore, universal precautions should be used in specimen handling (gloves, lab garments, avoid aerosol production, etc.). To determine HbA1c, a hemolyte must be prepared for each sample.

1. Dispense 1ml Hemolysis Reagent into tubes labeled: Control, Patients, etc. Note: Plastic or glass tubes of appropriate size are acceptable.

2. Place 20ul of well mixed whole blood into the appropriately labeled lysate reagent tube. Mix.

3. Allow to stand for 5 minutes or until complete lysis is evident. Hemolysates

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may be stored up to 10 days at 2-8°C.

Storage and Stability

1. All reagents are stable to the expiration date stated on the labels. Do not use the reagents past their expiration date.
2. R1 and R2 are stable for at least one month after opening stored at 2-8°C.
3. Hemoglobin A1c in whole blood collected with EDTA is stable for one week at 2-8°C.

Interferences

1. Bilirubin to 50mg/dL, ascorbic acid to 50mg/dL, triglycerides to 2000mg/dL, carbamylated Hb to 7.5mmol/L and acetylated Hb to 5.0mmol/L do not interfere in this assay.
2. It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin 6.78g.
3. It has been reported that elevated levels of HbF may lead to underestimation of HbA1C. Also, it has been reported that labile intermediates (Schiff base) are not detected and do not interfere with HbA1c determination by immunoassay 5.
4. It has been determined that Hemoglobin variants HbA2, HbC and HbS do not interfere with this method.
5. Other very rare variants of hemoglobin (e.g. HbE) have not been assessed.

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1. Pipettes to dispense 20 µl and 1 ml and Test Tubes to hold 1.02 ml.
2. Hemoglobin A1c calibrator set, control set, and Hemolysis rgt for 120 ml kit.

Procedure

(automated-Hitachi 717)
TEST NAME HbA1c
ASSAY CODE [1-POINT][50][0]
SAMPLE VOLUME [5] [3]
R1 VOLUME [180] [50] [NO]
R2 VOLUME [60] [20] [NO]
WAVELENGTH [] [660]
CALIBRATION [NONLINEAR] [4] [5]
STD (1) CONC-POS [0.0] [1]
TD (2) CONC-POS nm STD (3) CONC-POS no STD (4) CONC-POS nm STD
(5) CONC-POS [] [5]
STD (6) CONC-POS -
SD LIMIT [999]
DUPLICATE LIMIT [1000]
SENSITIVITY LIMIT [0]
ABS LIMIT (INC/DEC) [32000] [INCREASE]
PROZONE LIMIT [] []
EXPECTED VALUE [] []
PANIC VALUE [] []
INSTRUMENT FACTOR [1.0]
* Use Saline for the 0.0 Calibrator
** Input the values of the calibrator set being used TM Hitachi 717 is Ltd., Japan a
Registered Trademark of Nissei Sangyo Co.

Limitations

1. This assay should not be used for the diagnosis of diabetes mellitus.
2. Patient specimens should always be assayed using a calibration curve.
3. It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin 6-7.8g.
4. It has been reported that elevated levels of HbF may lead to underestimation of HbA1c and, that uremia does not interfere with HbA1c determination by immunoassay 10. It has been reported that labile intermediates (Schiff base) are not detected and therefore, do not interfere with HbA1c determination by immunoassay 5.
5. It has been determined that Hemoglobin variants HbA2, HbC and HbS do not interfere with this method.
6. Other very rare variants of hemoglobin (e.g. HbE) have not been assessed.

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TEST NAME	HbA1c
ASSAY CODE	[1-POINT][50][0]
SAMPLE VOLUME	[7] [3]
R1 VOLUME	[240] [50] [NO]
R2 VOLUME	[80] [20] [NO]
WAVELENGTH	[] [660]
CALIBRATION	[NONLINEAR] [4] [5]
STD (1) CONC-POS	[0.0] [1]
TD (2) CONC-POS	[**] [2]
STD (3) CONC-POS	[**] [3]
STD (4) CONC-POS	[**] [4]
STD (5) CONC-POS	[**] [5]
STD (6) CONC-POS	-
SD LIMIT	[999]
DUPLICATE LIMIT	[1000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[] []
EXPECTED VALUE	[] []
PANIC VALUE	[] []
INSTRUMENT FACTOR	[1.0]

* Use Saline for the 0.0 Calibrator

** Input the values of the calibrator set being used

TM Hitachi 717 is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan

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
Quality Control

The reliability of test results should be monitored whenever patient samples are assayed using a standard and quality control materials analyzed in the same manner employed for the unknowns. We suggest the use of commercially available Hemoglobin A1c controls with an assayed range. If controls do not fall into the assayed range patient values from that run should not be reported. The run should be repeated, making sure that all mixing and handling instructions are strictly followed.

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Calculations / Results

HbA1c results for the unknowns and controls are determined using the prepared calibration curve. An example curve is illustrated below:

HbA1c (%)

Expected Values¹¹

Recommended Values: less than 6% for a non-diabetic, less than 7% for glycemic control of a person with diabetes.

Each laboratory should establish its own expected values. In using Hemoglobin A1c to monitor diabetic patients, results should be interpreted individually. That is, the patient should be monitored against him or herself. There is a 3-4 week time lag before Hemoglobin A1c reflects changes in blood glucose level.

Performance

1. Linearity: The Hemoglobin A1c assay range is 2.0%-16.0%.

2. Comparison: A study using 40 human specimens between this Hemoglobin A1c procedure and an automated HPLC procedure (Tosoh) yielded a correlation coefficient of 0.988 and a linear regression equation of $y = 1.050x - 0.481$. ($S_{yx} = 0.332$)

3. Precision

Within Run: The within run precision was established by assaying two blood samples following NCCLS protocol EP5 on a Hitachi 917.

	Mean	S.D.	C.V.%
Low	5.48	0.078	1.4
High	10.28	0.176	3

Day to Day: The between day precision was established by assaying two blood samples following NCCLS protocol EP5 on a Hitachi 917.

	Mean	S.D.	C.V.%
Low	5.48	0.152	2.77
High	10.28	0.275	2.68

4. Sensitivity

Sensitivity was investigated by reading the change in absorbance at 660nm for a saline sample and a whole blood sample with a known concentration. Ten replicates of each sample were performed. The results of this investigation indicated that on the analyzer used (Hitachi 717), the HbA1c reagent

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showed little or no drift on the zero sample. Under the reaction conditions described, a 0.073 absorbance change is approximately equivalent to 1.0% HbA1c.

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INTENDED USE

For the *In Vitro* quantitative determination of glucose in human serum or plasma.

INTRODUCTION

Glucose is the major carbohydrate present in the peripheral blood. The oxidation of glucose is the major source of cellular energy in the body. Glucose determinations are run primarily to aid in the diagnosis and treatment of diabetes mellitus. Elevated glucose levels are mainly associated with insulinemia or insulin-induced hypoglycemia. A number of secondary factors also can contribute to elevated blood glucose levels. These include pancreatitis, pituitary or thyroid dysfunction, renal failure and liver disease. An enzymatic approach for glucose determination involves hexokinase coupled with glucose-6-phosphate dehydrogenase. A revision of this approach is proposed by the U.S. Center For Disease Control as the reference method for glucose and forms the basis of the reagent for glucose. 4

PRINCIPLE

Glucose + ATP $\xrightarrow{HK}$ G-6-P + ADP
G-6-P + NAD $\xrightarrow{G6PDH}$ 6-Phosphogluconate + NADH

The enzymatic hexokinase (HK) catalyzes the reaction between glucose and adenosine triphosphate (ATP) to form glucose-6-phosphate and adenosine diphosphate (ADP). In the presence of NAD, the enzyme glucose-6-phosphate dehydrogenase (G6PDH) oxidizes glucose-6-phosphate to 6-phosphogluconate. The increase in NADH concentration is directly proportional to the glucose concentration and can be measured spectrophotometrically at 340 nm.

REAGENT

Each kit contains four of 120 ml of Glucose Reagent (4 x 120 ml)

Instruction Insert.

Glucose standard (100 mg/dl)

REAGENT COMPOSITION

Adenosine triphosphate: 3.8 mmol/L

NAD⁺: 2.7 mmol/L

Hexokinase: 2000 IU/L

Glucose-6-phosphate dehydrogenase: 3000 IU/L

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PRECAUTIONS

1. This reagent is for *In Vitro* diagnostic use only.
2. The reagent contains sodium azide as a preservative. Do not ingest. Avoid skin and eye contact.
3. All specimens and controls should be handled as potentially infectious, use safe laboratory procedures (NCCLS M29-T2).

REAGENT PREPARATION

The reagent is ready to use.

REAGENT STORAGE

Reagent is stable until the expiration date on its respective label, when properly stored at 2 - 8°C and protected from light. Reagent should appear clear and colorless.

REAGENT DETERIORATION

Do not use if:

1. Reagent has an absorbance greater than 0.30 when measured against water at 340 nm.
2. The reagent fails to recover stated control values or meet stated linearity.
3. The reagent develops turbidity, indicating contamination.

SPECIMEN COLLECTION

1. Either serum or plasma may be used.
2. Plasma or serum samples without preservatives should be separated from the cells or clot within a half hour of being drawn.
3. Glucose in separated unhemolyzed serum is generally stable up to eight hours at 25°C and up to 72 hours at 4°C.
4. Glycolysis can be inhibited by collecting the specimen in sodium fluoride. Glucose in a sodium fluoride-oxalate mixture is reported to be stable up to 24 hours at 25°C.

INTERFERING SUBSTANCES

Grossly lipemic or icteric sera may cause falsely elevated glucose values, and may require the use of a serum blank. Young et al. give a complete list of drug and other substances that may affect glucose values. 6

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1. Accurate pipetting devices.
2. Timer.
3. Test tubes and rack.
4. Spectrophotometer capable of reading at 340 nm.

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MANUAL PROCEDURE

1. Appropriately label tubes: reagent blank, standard, patient sample, etc.
2. Pipette 1.0 ml of the glucose reagent into all tubes and prewarm to 37°C.
3. Add 0.005 ml (5 µl) of sample to respective tubes. Mix well. Incubate all tubes at 37°C for five minutes.
4. After incubation, zero the spectrophotometer with reagent blank at 340 nm.
5. Read and record the absorbance of all tubes.

* HTI-MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

LIMITATIONS

The procedure is linear to 500 mg/dl (27.8 mmol/L). Specimens above this limit must be diluted 1:1 with saline, re-run and the result multiplied by two (2) to compensate for the dilution.

CALIBRATION

An aqueous based standard is provided. The use of a serum-based calibrator is recommended.

QUALITY CONTROL

It is recommended that both normal and abnormal quality control sera be used routinely.

CALCULATIONS

Glucose results are expressed as mg/dl (mmol/L).

Abs. = absorbance at 340 nm

$$\frac{\text{Sample Abs.}}{\text{Standard Abs.}} \times \text{Conc. of standard} = \text{Conc. of glucose (mg/dl)}$$

$$\text{Example: Sample Abs.} = 0.155$$

$$\text{Standard Abs.} = 0.164$$

$$\frac{0.155}{0.164} \times 100 = 93 \text{ mg/dl}$$

Note: To convert the results into SI units (mmol/L), multiply the result (mg/dl) by 0.0556.

EXPECTED VALUES

Normal range is reported to be 65-100 mg/dl. This range should serve only as a guideline. It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations. In a study of 50 samples, the expected values were found to be 65-112 mg/dl.

PERFORMANCE CHARACTERISTICS

1. **Linearity:** 500 mg/dl.
2. **Sensitivity:** Based on an instrument resolution of 0.001 absorbance, this procedure has a sensitivity of 0.30 mg/dl.
3. **Comparison:** Glucose (HK) was compared to a commercially available glucose hexokinase method with the resulting linear regression equation of $Y = 0.999X + 0.44$. Coefficient of correlation $R^2 = 0.999$. Thirty-two patient sera and controls ranging from 46 mg/dl to 376 mg/dl were assayed by the two methods.
4. **Precision:**

	Within Run	
	Serum 1	Serum 2
Mean (mg/dl)	93	301
Std. Deviation (mg/dl)	1.58	2.59
C.V. (%)	1.70	0.90

	Run to Run	
	Serum 1	Serum 2
Mean (mg/dl)	93	300
Std. Deviation (mg/dl)	1.12	2.91
C.V. (%)	1.20	1.00

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REFERENCES

1. Cooper, G.R., *CRC Crit. Rev. Clin. Lab. Sci.* 4:101 (1973).
2. Henry, J.B., "Clinical Diagnosis and Management by Laboratory Method" W.B. Saunders and Company, Philadelphia, PA, p. 153 (1979).
3. Barthelmei, W., and Czok, R., *Klin. Wochenschr.*, 40:585 (1962).
4. A Proposed Method for Determining Glucose Using Hexokinase and Glucose-6-phosphate Dehydrogenase, Public Health Service, Center for Disease Control (1976).
5. Tietz, N.W., *Fundamentals of Clinical Chemistry*, 2nd. Ed., W.B. Saunders Co., Philadelphia, PA 243 (1976).
6. Young, D.S. et al., *Clin. Chem.* 21:5 (1975).

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htmed.com

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INTENDED USE

For the quantitative determination of total glucose in serum.


INTRODUCTION

Glucose is the major carbohydrate present in the peripheral blood. The oxidation of glucose is the major source of cellular energy in the body. Glucose determinations are run primarily to aid in the diagnosis and treatment of diabetes mellitus. Elevated levels glucose levels may be associated with pancreatitis, pituitary or thyroid dysfunction, renal failure and liver disease, whereas low glucose levels may be associated with insulinoma, hypopituitarism, neoplasms, or insulin induced hypoglycemia.

Early enzymatic methods for glucose determination involved glucose oxidase to catalyze the oxidation of glucose. Keston modified this method in the early 1950's using glucose oxidase/peroxidase enzyme system and o-dianisidine chromogen system.³ Since then, various alternative chromogen systems have been proposed. The trinder method replaces carcinogenic o-dianisidine with phenol plus 4-aminoantipyrine.⁴ This method is less influenced by interfering substances and does not suffer from the many drawbacks of earlier methods.

PRINCIPLE

The enzymatic reaction sequence employed in the assay of glucose is as follows:

β -D-Glucose is oxidized by glucose oxidase to produce D-gluconic acid and hydrogen peroxide. The hydrogen peroxide is then oxidatively coupled with 4-aminoantipyrine and phenol substitute, p-HBS, in the presence of peroxidase to yield a red quinonimine dye. The amount of colored complex formed is proportional to glucose concentration and can be photometrically measured.

REAGENT COMPOSITION

1. Glucose (Liquid) Reagent: Glucose Oxidase 15 IU/ml, Peroxidase (horseradish) 1.2 IU/ml, 4-Aminoantipyrine 0.2 mM, Phenol 4 mM, non-reactive ingredients and preservatives.
2. Glucose Standard: 100 mg/dl β -D-glucose in aqueous solution.

REAGENT PREPARATION

All reagents are ready for use.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use. **CAUTION:** *In vitro* diagnostic reagents may be hazardous. Handle in accordance with good laboratory procedures which dictate avoiding ingestion, and eye or skin contact.
2. Specimens should be considered infectious and handled appropriately.
3. Use distilled or deionized water where indicated.

REAGENT STORAGE AND STABILITY

Both liquid reagent and standard should be stored at 2 - 8°C. The reagent may be used until the expiration date indicated on the package label.

REAGENT DETERIORATION

The reagent should be discarded if:

1. Turbidity has occurred; turbidity may be a sign of contamination.
2. The reagent fails to meet linearity claims or fails to recover control values in the stated range.

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SPECIMEN COLLECTION

1. Test specimens should be serum and free from hemolysis.
2. Serum must be separated from the clot promptly since the rate of glucose decrease is approximately 7% per hour in whole blood.
3. Glucose in serum or plasma is stable for twenty-four (24) hours when stored at 2 - 8° C.

INTERFERING SUBSTANCES

Grossly lipemic or icteric sera will cause false glucose values and require the use of a serum blank. 5 Add 0.02 ml (20 µl) of patient sera to 3.0 ml distilled water and read against a water blank. Subtract this absorbance from the patient test absorbance to correct for the lipemia or icterus. Young, et al. give a comprehensive review of drug interferences.

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1. Pipettes to accurately measure required volumes
2. Test tubes/rack
3. Timer
4. 37 °C heating block or water bath
5. Spectrophotometer capable of accurately measuring absorbances at 500 nm

GENERAL INSTRUCTIONS

The reagent for Glucose Liquid Reagent is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

AUTOMATED PROCEDURE

Please refer to appropriate application manual available.

MANUAL PROCEDURE

1. Label tubes: blank, standard, control, patient, etc.
2. Pipette 1.5 ml of working reagent to all tubes and place in a 37°C heating bath for at least five (5) minutes.
3. Add 0.01 ml (10 µl) of sample to respective tubes, mix and incubate at 37°C for exactly ten (10) minutes.
4. After incubation, zero spectrophotometer with the reagent blank. Read and record the absorbances of all tubes at 500 nm (Wavelength range: 500 - 520 nm).

* HTI MULTI-PURPOSE CALIBRATOR MAY BE USED TO

REPLACE STANDARD

NOTE: If the spectrophotometer being used requires a final volume greater than 1.5 ml for accurate reading, use 0.02 ml (20 µl) of sample to 3.0 ml of reagent. Perform the test as described above.

CALIBRATION

The procedures are calibrated with the standard solution which is included with each series of tests. Its absorbance is used to calculate results. It is recommended to establish a linearity curve up to 500 mg/dl with other available commercial standard solutions to verify the performance of the instruments and reagents.

LIMITATIONS

The reagent is linear to 500 mg/dl glucose. Samples that have glucose values greater than 500 mg/dl should be diluted with water 1:1, re-assayed and the results multiplied by

CALCULATIONS

(A = Absorbance)

$$\frac{A(\text{patient})}{A(\text{standard})} \times \text{Concentration of standard (mg/dl)} = \text{Concentration of unknown (mg/dl)}$$

Example:

A (patient) = 0.37, A (standard) = 0.28
Concentration of standard = 100 mg/dl

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HTI-G520-1000

$$\frac{0,37}{0,28} \times 100 = 132 \text{ mg/dl}$$

SI UNITS: To obtain results in SI units (mmol/L), multiply your result in mg/dl by ten (10) to convert dl to liter and divide the value by 180, the molecular weight of glucose.
 $\text{mg/dl} \times 10 = \text{mg/dl} \times 0,0556$
 180

Example: $132 \text{ mg/dl} \times 0,0556 = 7,34 \text{ mmol/L}$

QUALITY CONTROL

It is recommended that high and low values of glucose controls be included in each set of assays. Commercially available control material with established glucose values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction or procedural errors.

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PERFORMANCE CHARACTERISTICS

- Linearity:** 500 mg/dl
- Sensitivity:** An absorbance change of 0.001 at 500 nm corresponds to 0.5 mg/dl under the stated condition of this assay system.
- Comparison:** A comparison between this reagent in liquid form and the powder form produced a regression equation of $y = 1,00x + 2,54$ (N=64) with a coefficient of correlation of 0.99.

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Intended Use

For the quantitative, kinetic determination of glucose-6-phosphate dehydrogenase (G6PD) in blood at 340 nm. For *in vitro* diagnostic use only. Clinical Significance: G6PD assays are most commonly performed to determine deficiency of G6PD, which is widely prevalent throughout the world. It has been determined that G6PD deficiency in red cells is the basis for certain drug-induced hemolytic anemias. This type of susceptibility to drug-induced hemolysis is often called "primaquine sensitivity" because studies which led to its characterization were made during investigations of the hemolytic properties of this antimalarial compound.

Summary

Glucose-6-phosphate dehydrogenase (G6PD, D-glucose-6-phosphate oxidoreductase, EC 1.1.1.49) catalyzes the first step in the pentose phosphate shunt, oxidizing glucose-6-phosphate (G-6-P) to 6-phosphogluconate (6-PG) and reducing NADP to NADPH. This procedure is a modification of the spectrophotometric methods of Kornberg and Horecker² and of Lohr and Waller³, involving the following reaction:

G-6-P + NADP⁺ → 6-PG + NADPH + H⁺
Nicotinamide adenine dinucleotide phosphate (NADP) is reduced by G6PD in the presence of G-6-P. The rate of formation of NADPH is proportional to the G6PD activity and is measured spectrophotometrically as an increase in absorbance at 340nm. Production of a second molar equivalent of NADPH by erythrocyte glucose-6-phosphate dehydrogenase (6-PGDH) according to the reaction:

6-PG + NADP⁺ → Ribulose-5-Phosphate + NADPH + H⁺ + CO₂
is prevented by use of maleimide, an inhibitor of 6-PGDH.

Reagents

G6PD R1 Reagent: Reconstituted reagent will contain NADP, 1.5 mM, and maleimide, 12 mM. Also contains buffer, stabilizer and lysing agent.
G6PD R2 Reagent: Glucose-6-Phosphate, 1.05 mM, buffer and magnesium salt. Sodium azide added as preservative.

Precautions

1. These reagents are for *in vitro* diagnostic use only.
2. Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state and federal laws.
3. R1 Reagent is HARMFUL. May cause sensitization by inhalation and skin contact. Wear suitable protective clothing.
4. R2 Reagent contains sodium azide which may react with lead and copper plumbing to form highly explosive metal azides. Avoid azide accumulation.

Reagent Preparation

1. R1 reagent is prepared by reconstituting with the volume of deionized water indicated on label or application sheet. Swirl gently and invert several times to dissolve contents. Wait 2-3 minutes and mix again. NOTE: For manual use, see reagent preparation instructions listed in MANUAL PROCEDURE section.
2. R2 reagent is supplied ready to use.

Storage and Stability

1. When unopened R1 reagent vials and the R2 reagent are stored at 2-8°C. They are stable until the expiration date on the labels.
2. Reconstituted R1 reagent is stable for 8 hours at room temperature (18-26°C) or 5 days refrigerated (2-8°C).

Optional Reagents

G6PD controls: Lyophilized controls with G6PD in a stabilized human red cell hemolysate base. High Technology catalog number G7583-CTL.
G6PD lysing reagent: Triton X-100, 0.05% v/v. For use with discrete analyzer applications. High Technology catalog number G7583-LYS.

Specimen Collection and Storage

1. It is recommended that specimen collection be carried out in accordance with NCCLS document M29-T2.
2. Whole blood collected in EDTA, heparin or acid-citrate-dextrose (ACD) is

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satisfactory 4-8

3. Red cell G6PD is stable in whole blood for one week refrigerated (2-8°C), but is unstable in red cell hemolysates. 9
4. Freezing of blood is not recommended. 4
5. Since activity is reported in terms of grams hemoglobin or the number of red blood cells, the hemoglobin concentration or red cell count must be determined prior to performing the G6PD assay. The integrity of erythrocytes collected in ACD is preserved even after prolonged storage so that obtaining accurate red cell counts usually poses no problem. 6 However, red cell counts on specimens collected in heparin become unreliable after about 2 days. 6 Thus, for heparinized samples, results are best reported in terms of hemoglobin concentration.

Interfering Substances

1. Copper completely inhibits G6PD at a concentration of 100 $\mu\text{mol/L}$, and sulfate ions (0.005 mol/L) decrease observed levels of G6PD activity. 10
 2. Certain drugs and other substances are known to influence circulating levels of G6PD. 11
 3. Reticulocytes have higher G6PD levels than mature red cells. It is recommended that assays not be performed after a severe hemolytic crisis, since G6PD levels may appear falsely elevated. Under those conditions, detection of deficiency may require family studies. Testing may be performed after the level of mature red cells has returned to normal.
 4. Under normal circumstances, activity contributed by leukocytes, platelets and serum is relatively small. However, in cases of extreme anemia, grossly elevated white counts, or very low levels of red cell G6PD activity, the contribution to the total made under these conditions may be significant. See "Use of Buffy-Coat-Free Samples" section.
- Automated Analyzer Applications
Application procedures are available for various automated instruments. Please contact High Technology Inc. for more information.

Materials Provided

See "reagents" section.

Materials Required But Not Provided

1. Spectrophotometer capable of measuring at 340 nm with temperature controlled cuvette compartment (water bath or incubator may be used instead).
2. Pipetting devices for delivery of volumes required for the assay. Cuvettes with optical properties suitable for use at 340 nm.
3. Equipment and reagents for determining hemoglobin concentration or performing a red cell count. High Technology Inc.

Manual Procedure

Prepare working R1 reagent by adding lyse as the diluent instead of DH2O. Add volume stated on the R1 vial. This reagent can now be used as stated below. **NOTE: Do not use DH2O to reconstitute the R1 vial for manual procedure.**

The temperature of the reaction must be maintained at 37°C or some other constant temperature (see "Temperature Correction" section).

1. Prepare reaction mixture:
 - a. To a labeled cuvette, add 1.0 ml R1 Reagent.
 - b. Add 0.01 ml blood and mix thoroughly to completely suspend erythrocytes. Let stand at room temperature (18-26°C) for 5-10 minutes.
 - c. Add 2.0 ml R2 Reagent and mix gently by inverting several times.

Proceed to step 2.

2. Place cuvette in constant temperature cuvette compartment or water bath and incubate for approximately 5 minutes.
3. Read and record absorbance (A1) of TEST at 340 nm vs water. (If using a water bath or incubator, return cuvette to it.)
4. Exactly 5 minutes later, read and record absorbance (A2).

5. To determine G6PD activity, refer to "Calculations" section.

Calibration

The procedure is standardized on the basis of the millimolar absorptivity of NADPH, which is 6.22 at 340 nm. Measurement of the rate of increase in absorbance (A2) at 340 nm serves to quantitate enzymatic activity.

Quality Control

Reliability of test results should be monitored by use of control materials with known values within each run. High Technology Inc. glucose-6-phosphatase

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dehydrogenase controls are suitable for this purpose. (Catalog number

G7583-CTL) It is recommended that each laboratory establish its own frequency of control determination.

Calculations

ΔA per min = $A_2 - A_{15}$
G6PD activity can be expressed as either U/g hemoglobin (Hb) or as U/1012 erythrocytes (RBC).

$$\text{G6PD (U/g Hb)} = \Delta A \text{ per min} \times \frac{100 \times 3.01}{0.01 \times 6.22 \times \text{Hb (g/dl)}} \times \text{TCF}$$

$$= \Delta A \text{ per min} \times \frac{4839}{\text{Hb (g/dl)}} \times \text{TCF}$$

Where: 100 = Factor to convert activity to 100 ml
3.01 = Total reaction volume (ml)
0.01 = Sample volume (ml)
6.22 = Millimolar absorptivity of NADPH at 340 nm
Hb (g/dl) = Hemoglobin concentration for each specimen
TCF = Temperature Correction Factor (1 at 37°C)
or G6PD (U/1012RBC) = $\Delta A \text{ per min} \times \frac{3.01 \times 10^{12} \times \text{TCF}}{0.01 \times 6.22 \times (N \times 10^6) \times 1000}$

Where: 3.01 = Total reaction volume (ml)
1012 = Factor for expressing activity in 1012 cells
0.01 = Sample volume (ml)
6.22 = Millimolar absorptivity of NADPH at 340 nm
 $N \times 10^6$ = Red cell count (red cells/mm³) determined for each specimen
1000 = Conversion of red cell count from mm³ to ml
TCF = Temperature correction factor (1 at 37°C)

This equation reduces to:
G6PD (U/1012 RBC) = $\Delta A \text{ per min} \times \frac{48,390}{N} \times \text{TCF}$

Where: N = red cell count divided by 10⁶
*TCF = temperature correction factor (1 at 37°C)
* Additional temperature correction factors are available upon request.

Example

Assay of a specimen which had a red cell count of 4.6 x 10⁶/mm³ and a hemoglobin concentration 15.2 g/dl resulted in a ΔA per min at 37°C of 0.028.
G6PD (U/g Hb) = $0.028 \times \frac{4839}{15.2} = 8.9$

$$\text{G6PD (U/1012 RBC)} = 0.028 \times \frac{48,390}{4.6} = 295$$

NOTE: If ΔA per min is greater than 0.060, repeat determination using 5 μ l blood and multiply results by 2.

Use of Buffy-Coat-Free Sample

Under normal circumstances G6PD activity contributed by leukocytes, platelets and serum is relatively small. However, as reported by Echler¹² and others¹³, more accurate measurement of red cell G6PD activity, especially in the presence of anemia and/or leukocytosis, can be achieved by using buffy coat-free blood samples for assay. Thus, in case of a borderline value obtained with whole blood, it may be preferable to repeat the assay on a buffy coat-free sample.

Temperature Correction

When temperature is 37°C, no temperature correction factor (TCF) is required in the calculations. If assay is performed at another temperature a TCF must be used.¹⁴

Cuvette Temperature	TCF
25°C	1.98
30°C	1.37

Unit Definition

One International Unit (U) is that amount of G6PD activity that will convert 1 micromole of substrate per minute under the conditions specified in this insert.

Expected Values¹

A recommended reference range for G6PD measured at 37°C is:

$$12.1 \pm 2.09 \text{ U/g Hb}$$

$$351 \pm 60.6 \text{ U/1012 RBC}$$

Values for newborns may range somewhat higher. It is highly recommended that each laboratory establish its own expected range.

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
Performance Characteristics

Assay Range The maximum G6PD activity which may be measured by this procedure is approximately 21.0 U/g Hb or 609 U/1012 RBC.

Theoretical

Observed Data	Theoretical Result	Percent Recovery
2.78 U/g Hb	2.78 U/g Hb	100.0%
5.29 U/g Hb	5.56 U/g Hb	95.1%
10.80 U/g Hb	11.12 U/g Hb	97.1%
20.69 U/g Hb	22.24 U/g Hb	93.0%

Precision:

Precision studies were performed on a Roche Cobas Mira following the guidelines contained in NCCLS document EPS-T2.15. The data is presented in units that an automated analyzer will produce for G6PD activity (U/L). It is highly recommended that precision of the assay be verified on each analyzer before use.

Within Day (n=20)

Mean	S.D.	C.V.
257	23.7	9.2%
658	18.3	2.8%
1939	48.0	2.5%

Day to Day (n=20)

Mean	S.D.	C.V.
269	30.8	11.4%
700	28.7	4.1%
2014	43.0	2.1%

Sensitivity:

Assuming the limit of sensitivity to represent a change in absorbance at 340nm of 0.001 per minute, a G6PD activity of 0.4 U/g Hb or 11 U/1012 RBC may be detected using this procedure (assuming a hemoglobin concentration of 12.0 g/dL and a red cell count of $4.5 \times 10^6/\text{mm}^3$).

Specificity:

The oxidation of glucose-6-phosphate by G6PD is specific. Any non-specific formation of NADPH due to oxidation of other substrates by endogenous enzymes occurs during the preincubation period. 6-Phosphogluconate dehydrogenase is completely inhibited by maleimide in the reagent system.

Correlation:


A comparison study between the High Technology, Inc. method and that of Sigma Diagnostics yielded a linear regression equation with $y = 0.97x + 0.07$ and a correlation coefficient of 0.994.

References

- Burtis, C.A., Ashwood, E.R., Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, pp. 1645-1650, 1988.
- Kornberg, A., Horecker, B.L. Glucose-6-Phosphate Dehydrogenase. IN Methods in Enzymology. S.P. Colowick, N.O. Kaplan, Editors, Vol. 1, Academic Press, New York, p 323, 1955.
- Lohr, G.W., Waller, H.D.: Glucose-6-phosphate Dehydrogenase. IN Methods of Enzymatic Analysis. H.U. Bergmeyer, Editor, Academic press, New York, p. 636, 1974.
- Kachmar, J.F., Moss, D.W. Enzymes. IN Fundamentals of Clinical Chemistry. N.W. Tietz, Editor, Saunders, Philadelphia, pp. 666-672, 1976.

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5. WHO Technical Report Series No. 366, Standardization of Procedures for the study of Glucose-6-Phosphate Dehydrogenase, 1967.
6. Lowe, M. L., Stella, A. F., Mosher, B. S., Gh, J. B., Demetriou, J. A. Microfluorimetry of Glucose-6-phosphate dehydrogenase and 6-phosphogluconate Dehydrogenase in red cells. Clin Chem 16:440, 1972.
7. Bishop, C., Assay of glucose-6 phosphate dehydrogenase (EC 1.1.1.49) and Glucose-6-phosphate dehydrogenase (EC 1.1.1.44) in red cells. J Lab Clin Med 68:149, 1966.
8. Beutler, E., Blume, K. G., Kaplan, C., Lohr, W., Ramont, S., Valentine, W. N.: International committee for standardization in haematology: Recommended screening test for glucose-6-phosphate dehydrogenase (G6PD). Br J of Haem, 43:469-477, 1979.
9. Siene, E. A. Red Cell Enzyme Deficiencies: A Review. Am J Med Tech 38:454, 1972.
10. Boulard M, Blume KG, Beutler E. The effects of copper on red cell enzyme activities. J. Clin Invest, 51, 459 (1972).
11. Young, D. S., Pestaner, L. C., Gibberman, V. Effects of drugs on clinical laboratory tests. Clin Chem 21: 302D, 1975.
12. Echler, G. Determination of glucose-6-phosphate dehydrogenase. Am J Med Technol 49:259, 1983.
13. Morelli, A., Benatti, U., Lenzerini, L., Sparatore, B. et al. The interference of leukocytes and platelets with measurement of glucose-6-phosphate dehydrogenase activity of erythrocytes with low activity variants of the enzyme. Blood 58: 642, 1981.
14. Beutler, E., et al. International Committee for Standardization in Haematology: Recommended methods for red-cell enzyme analysis. Br J Haematol., 35:331-340, 1977.
15. NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Edition, 1992.

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Homocysteine Reagent Set Гомоцистеин

Homocysteine

Type: reagent
Method: Enzyme Cycling/
Format: Liquid
Linearity: 3-50 $\mu\text{mol/l}$
Sample: serum

Kit composition:
Reagent P1: 1x54 ml
Reagent P2: 1x11 ml
Reagent P3: 1x8 ml
Kt (190 Tests)

Product Features:

The Diazyme enzymatic homocysteine assay reagent is a liquid stable system, which provides a fast, accurate cost effective alternative to cumbersome and expensive immunochemical methods. A choice of convenient instrument specific packaging options eliminate reagent transfer for Roche Hitachi, Dade Dimension®, Dade Expand®, and Beckman Synchron CX, LX and DxC analyzers. Accurate and reliable, the assay features excellent correlation with HPLC and Immunochemical methods with a linear range of 3-50 $\mu\text{mol/L}$, and intra and inter %CV values of < 5%. The assay has no significant interference from lipemic and hemolytic substances and eliminates the "carry-over" interference with iron or lipase reagents that can occur with other methods. More significantly, the assay is not interfered by the significantly, the assay is not interfered by the Enzymatic Homocysteine (3 Reagent)

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
2-10°C.

Intended use

The DCP Immunological Test System is an in vitro device that consists of reagents used to quantitatively measure, by immunochemical techniques, des-gamma-carboxy-prothrombin (DCP) in human serum. The test is intended for in vitro diagnostic use as an aid in the risk assessment of patients with chronic liver disease for development of hepatocellular carcinoma (HCC) in conjunction with other laboratory findings, imaging studies and clinical assessment.

Summary and explanation of the test

Prothrombin is a vitamin K dependent blood coagulation factor that is formed in the liver. It contains carboxy-glutamic acid residues on its amino-terminal domain, which are synthesized from glutamic acid residues by vitamin K dependent glutamyl carboxylase in a post-translational process. In the case of deficiency of vitamin K or the ingestion of vitamin K antagonists DCP is found in patients.

Principle of the method

The DCP assay is an easy-to-use system with all reagents in a single cartridge and each assay performed on a single, disposable "chip" using microfluidic electrophoretic separation. After placing sample, reagent cartridge, Wash Solution and Chip Cassette in the instrument, the buffers, antibody solutions and sample are automatically dispensed into appropriate chip wells. All reactions, separations and detection occur on a microfluidic chip.

Reagents

The DCP kit contains a reagent cartridge and an adapter. The reagent cartridge contains the following buffers and antibodies for conducting 100 tests:

- (1) Electrophoresis Buffer (R1) 5.4 mL 100 mmol/L Tris buffer, pH 7.5
- (2) Electrophoresis Buffer (R2) 4.4 mL 40 mmol/L Tris buffer, pH 8.0
- (3) Electrophoresis Buffer (R3) 2.6 mL 75 mmol/L Tris buffer, pH 7.1
- (4) Labeled Antibody Solution (C1) 0.77 mL 75 mmol/L Good's buffer, pH 6.0
- (5) 200 nmol/L Anion-conjugated anti human DCP antibody (mouse monoclonal antibody) (DNA-Fab' (DCP))
- (6) Labeled Antibody Solution (C2) 0.92 mL
- (7) 50 mmol/L Phosphate buffer, pH 5.55
- (8) 1000 nmol/L Fluorescent dye labeled anti human prothrombin antibody (mouse monoclonal antibody) (Dye-Fab' (prothrombin))
- (9) Fluorescent Dye Solution (FD) 1.4 mL
- (10) 50 mmol/L Good's buffer, pH 6.0

Store the reagent cartridge at 2-10°C (Do not freeze)

Specimen collection and preparation

- (1) Use serum as a specimen.
- (2) If immediate analysis is not possible, store specimen at -80°C. DCP concentration in serum is stable for 4 years frozen at -80°C. DCP in serum is stable for 3 weeks at -20°C. DCP in serum is stable for 1 week at -4°C.

Reagent preparation

Reagents: Use as supplied. Unopened reagents are stable until expiration date printed on the label.

Test procedure

Dispense reagents and sample to chip

Reaction with Dye-Fab'

The final results are automatically calculated and printed out or sent out to the host computer. The results are reported as a concentration of DCP.

Calibration

Calibration is required when reagent cartridge is opened. Calibration curve is automatically produced by plotting fluorescence intensity of the immunocomplex peak area vs. DCP concentrations of the calibrators. The calibration curve is stable for 30 days. For a detailed description on the calibration, refer to Section 4.2 "Calibration Operation" of the Instruction Manual.

Quality control

A quality control program is recommended for all clinical laboratories. Daily analysis using DCP Controls L and H is recommended for monitoring the performance of the procedure. The values obtained for the controls should fall within 15% for DCP across upper reportable range (> 1 ng/mL) from the assigned values or within 20% for DCP across lower reportable range (< 1 ng/mL) from the assigned values.

Reportable range

The reportable range of DCP concentration is 0.1 - 950 ng/mL.

Expected values

Less than 7.5 ng/mL.

Precautions

For in vitro diagnostic use only.

Handle serum and blood specimens as potentially infectious samples and follow

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LINERITY: 0.5 - 30 µg/mL (or value of highest calibration point)

SAMPLE: Serum, Plasma

KIT COMPOSITION:

R1: 1 X 16 mL buffer

R2: 1 X 8.5 mL antibody

Tests / Kit: 80 (Hitachi 717)

(Can Vary Based On Analyzer Model)

METHOD

In this method uses polystyrene microparticles coated by covalently bound monoclonal anti-D-dimer antibodies. If the sample contains D-dimer, the formation of complexes with these antibodies takes place. The formation of complexes leads to increase in absorbance of the reaction mixture, which is proportional to the concentration of antigen.

INTENDED USE:

The D-Dimer is an *in vitro* diagnostic reagent for the quantitative determination of cross-linked fibrin degradation products containing D-Dimer in human plasma or serum by immunoturbidimetric assay.

FOR *IN VITRO* DIAGNOSTIC USE.

REAGENT PREPARATION:

The reagents are ready to use. They do not need to be reconstituted. Mix before using by gently inverting the bottles. After opening, gently invert Reagent 2 once a week.

STORAGE AND HANDLING:

All reagents should be stored refrigerated (2-8°C) and protected from light. Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for one year from the date of manufacture, as indicated by the expiration date on the package and bottle labels. Opened reagents can be used for one month if stored at 2-8°C.

PERFORMANCE:

Precision

Precision - Within Run

Sample I	Serum Sample II	Sample III
N = 5	N = 5	N = 5
Mean = 2.29	Mean = 9.60	Mean = 22.04

S.D. = 0.06	S.D. = 0.14	S.D. = 0.36
C.V. = 2.7%	C.V. = 1.4%	C.V. = 1.6%

Sample I	Plasma Sample II	Sample III
N = 10	N = 10	N = 10
Mean = 3.55	Mean = 11.34	Mean = 21.12
S.D. = 0.09	S.D. = 0.12	S.D. = 0.12
C.V. = 2.5%	C.V. = 1.1%	C.V. = 0.6%

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TEST NAME	[FE]
ASSAY CODE	[2-POINT];[24];[50]
SAMPLE VOLUME	[15] [10]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [100] [NO]
WAVELENGTH	[700] [570]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[1] [1]
STD (3) CONC-POS	
STD (4) CONC-POS	

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STD (6) CONC-POS
SD LIMIT [0.1]
DUPLICATE LIMIT [200]
SENSITIVITY LIMIT [0]
ABS LIMIT (INC/DEC) [0] [INCREASE]
PROZONE LIMIT [0] [LOWER]
EXPECTED VALUE ["] ["]
PANIC VALUE ["] ["]
INSTRUMENT FACTOR [1.0]

* Indicates user defined parameter.

** Input appropriate bottle size.

Set K-Factor to 5160

Calculations

A = Absorbance

Std = Standard

$A_2 \text{ Test} - A_1 \text{ Test} \times \text{Conc.} = \text{Total Iron (ug/dl)}$

$A_2 \text{ Std} - A_1 \text{ Std} \text{ of Std}$

Example: $A_1 \text{ Test} = 0.08$

$A_2 \text{ Test} = 0.15$

$A_1 \text{ Std} = 0.00$

$A_2 \text{ Std} = 0.40$

Then: $\frac{0.15 - 0.08}{0.40 - 0.00} = \frac{0.07}{0.40} \times 500 = 0.175 \times 500 = 87.5 \text{ ug/dl}$

0.40 - 0.00 0.40

Calibration

Aqueous standards can be used to calibrate the procedure or an appropriate serum calibrator. The procedure should be calibrated according to the instrument manufacturer's instructions. If control results are found to be out of range, the procedure should be re-calibrated.

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Expected Values

Iron, Total = 60 - 150 ug/dl

It is strongly recommended that each laboratory determine the normal range for its particular population.

Performance

STD (5) CONC-POS
SD LIMIT [0.1]
DUPLICATE LIMIT [200]
SENSITIVITY LIMIT [0]
ABS LIMIT (INC/DEC) [0] [INCREASE]
PROZONE LIMIT [0] [LOWER]
EXPECTED VALUE ["] ["]
PANIC VALUE ["] ["]
INSTRUMENT FACTOR [1.0]

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HTI-1592-120

HTI-1592- 500

- | | | | | | | |
|--|-----|-----|-----|-----|-----|-----|
| 1. Linearity: 500 ug/dl | 110 | 3.0 | 2.7 | 109 | 5.1 | 4.7 |
| Samples with values above 500 ug/dl must be diluted 1:1 with normal saline, re-assayed and result multiplied by two. | 184 | 3.1 | 1.7 | 177 | 3.8 | 2.1 |
| 2. Comparison: A study performed between this procedure and a similar Serum Iron procedure resulted in a coefficient of correlation of .983 with a regression equation of $y = 0.79x + 0.13$ | 251 | 3.3 | 1.3 | 249 | 4.5 | 1.8 |
| 3. Precision: Total Iron | | | | | | |

Within Run			Run to Run		
Mean	S.D.	C.V.%	Mean	S.D.	C.V.%
110	3.0	2.7	109	5.1	4.7
184	3.1	1.7	177	3.8	2.1
251	3.3	1.3	249	4.5	1.8

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htimed.com

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Intended Use

For the quantitative determination of iron in serum. For *in vitro* diagnostic use only.

Method History

Iron exists in serum complexed with transferrin, a transport protein. Most early procedures for iron determination involved dissociation of the iron from the iron-protein complex, precipitation of the proteins, and then measurement of the iron content of the protein free filtrate.

Many chromagens have been used in the determination including thiocyanate *o*-phenanthroline, bathophenanthroline and TPTZ. In 1971, Persijn et al.¹ presented a method using the chromagen ferrozine, described by Stookey.² This method did not require protein precipitation and was more sensitive than previous methods. The present procedure is a modification of the Persijn method.

Principle

Serum Iron: Transferrin-bound iron is released at an acid pH and reduced from ferric to ferrous ions. These ions react with ferrozine to form a violet colored complex which is measured spectrophotometrically at 560nm. The absorbance measured at this wavelength is proportional to serum iron concentration.

Clinical Significance³

In most cases, both serum iron and TIBC values are necessary for greatest diagnostic significance. Low serum iron values are seen in chronic blood loss, insufficient intake or absorption of iron, and increased demand on the body stores (e.g. pregnancy). Elevated serum iron values are seen in increased red cell destruction, decreased red cell synthesis, increased iron intake, or increased iron stores release.

Increase in the TIBC may be due to increased production of apotransferrin (e.g. chronic iron deficiency) or an increased release of ferritin, as in hepatocellular necrosis.

Decreases in the TIBC can occur with cirrhosis and hemochromatosis due to a deficiency in ferritin, or in nephrosis due to loss of apotransferrin.

Reagents

1. IRON BUFFER REAGENT: Hydroxylamine hydrochloride 220mM in acetate buffer, pH 4.5 with surfactant.
2. IRON COLOR REAGENT: Ferrozine 16.7mM in hydroxylamine hydrochloride.

Precautions

1. All reagents are toxic. Do not pipette by mouth. Avoid all contact.
2. This reagent is for *in vitro* diagnostic use only.

Reagent Storage

Store all reagents refrigerated at 2-8°C.

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Reagent Deterioration

All reagents should be clear. Turbidity may indicate contamination and the reagent should not be used.

Specimen Collection and Storage

1. Fresh, unhemolyzed serum is the specimen of choice.
2. Serum should be separated as soon as clot has formed.
3. Heparinized plasma may be used but other anticoagulants should not be used to avoid possible iron contamination.⁴
4. Serum iron is reported to be stable for four days at room temperature (15-30°C) and seven days at 2-8°C.⁴

Interferences

1. Certain drugs and other substances are known to influence circulating iron levels. See Young, et al.⁵
2. Iron contained in hemoglobin does not react in this method, therefore, slight hemolysis will not interfere. However, gross hemolysis (pink or red specimens) will contribute to the absorbance measured at the wavelength used and should be avoided.³
3. To make tubes, pipettes, etc. iron free, they must be washed with hot, dilute (1:2) hydrochloric or nitric acid, followed by several rinsings with iron-free deionized or distilled water.

Materials Provided

1. Iron Buffer Reagent.
2. Iron Color Reagent.
3. Iron Standard

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1. Accurate pipetting devices
2. Test tubes/rack
3. Timer
4. Spectrophotometer able to read at 560 nm.
5. Iron-free deionized water.
6. Heating bath/block (37°C).

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Procedure (Manual)

Serum Iron

1. Label test tubes/cuvettes, "Blank", "Standard", "Control", "Sample", etc.
2. Add 2.5ml Iron Buffer reagent to all tubes.
3. Add 0.5ml (500ul) sample to respective tubes. Mix. Note: Add 500ul iron-free water to blank.
4. Zero spectrophotometer at 560nm with the reagent blank.
5. Read and record absorbances of all tubes. (A₁ reading).
6. Add 0.05ml (50ul) Iron color reagent to all tubes. Mix.
7. Place all tubes in heating bath at 37°C for 10 minutes.
8. Zero instrument at 560nm with reagent blank.
9. Read and record absorbances of all tubes. (A₂ reading).

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Calibration

The procedure is calibrated with an aqueous iron standard (500ug/dl) or an appropriate serum calibrator.

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Performance

- Linearity: 500 ug/dl
Samples with values above 500 ug/dl must be diluted 1:1 with normal saline, re-assayed and result multiplied by two.
- Comparison: A study performed between this procedure and a similar Serum Iron procedure resulted in a coefficient of correlation of .993 with a regression equation of $y = 1.02x + 7.0$
- Precision: Total Iron

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- Saunders, pp. 923-929, (1976).
4. Weissman, N., Pileggi, V.J., in Clinical Chemistry: Principles and Technics, 2nd Ed., R.J. Henry et al, editors, Hagerstown (MD), Harper & Row, pp. 692-693, (1974).
5. Young, D.S. et al, Clin. Chem. 21:1D, (1975).
- Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods, Philadelphia, W.B. Saunders, p. 1434, (1984).

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htimed.com

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KIT COMPOSITION

R1: 4 X 20 mL buffer
 R2: 2 X 10 mL antibody
 Tests / Kit: 285 (Hitachi 717)
 (Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human IgA, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm.

INTENDED USE

For the quantitative determination of human IgA in serum by immunoturbidimetric assay. Measurement of IgA aids in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents.

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KIT COMPOSITION

R1: 4 X 20 mL buffer
 R2: 4 X 20 mL antibody
 Tests / Kit: 320 (Hitachi 717)
 (Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human IgG, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm.

INTENDED USE

For the quantitative determination of human IgG in serum by immunoturbidimetric assay. Measurement of IgG aids in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents.
 FOR IN VITRO DIAGNOSTIC USE.

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KIT COMPOSITION

R1: 4 X 20 mL buffer
R2: 2 X 10 mL antibody
Tests / Kit: 285 (Hitachi 717)
(Can Vary Based On Analyzer Model)

METHOD

The method is based on the reaction of specific antibodies with samples containing human IgM, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm.

INTENDED USE

For the quantitative determination of human IgM in serum by immunoturbidimetric assay. Measurement of IgM aids in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents.
FOR IN VITRO DIAGNOSTIC USE.

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KIT COMPOSITION:

Option 1

R1: 1 X 13.5 mL buffer
R2: 1 X 5 mL antibody
Tests / Kit: 100 (Hitachi 917)
(Can Vary Based On Analyzer Model)

Option 2

R1: 2 X 13.5 mL buffer
R2: 2 X 5 mL antibody

Reagents (Liquid Stable)

R1: Buffer Reagent
Tris(hydroxymethyl)aminomethane, 100 mM, pH 8.2
R2: Latex Suspension
Anti-human insulin mouse monoclonal antibody, ~1 mg/mL

METHOD

The method is based on the reaction of specific antibodies with samples containing human insulin, which resulted in formation of the insoluble complex. The presence of the complex can be measured by optical density at 340 nm. The concentration of insulin in the sample can be determined using standard curve constructed by measuring the optical density of standards.

INTENDED USE:

For the quantitative determination of human insulin in serum and plasma by immunoturbidimetric assay
FOR IN VITRO DIAGNOSTIC USE

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Assay Procedure

An example of automated application:
Sample 12 μ L

↓
□ ← R1 (Buffer Reagent) 135 μ L

↓ 37°C, 5 min.

□ ← R2 (Latex Reagent) 50 μ L

↓ 37°C, 5 min.

Rate, 600 nm (main) / 800 nm (sub)

Note: Allow all reagents and specimens to warm to room temperature. Mix all reagents gently before using.

PERFORMANCE:

Precision

The within-run, between-run, and total precision for insulin assay was determined using packaged reagents, human serum samples, and a Roche / Hitachi 917 analyzer.

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Insulin Reagent Set Инсулин

Precision - Between Runs

Sample I	Sample II	Sample III
N = 10	N = 10	N = 10
Mean = 12.21	Mean = 25.9	Mean = 68.74
S.D. = 0.530	S.D. = 0.777	S.D. = 1.693
C.V. = 4.34%	C.V. = 3.00%	C.V. = 2.46%

Accuracy / Correlation

A comparison of Insulin assay and another company's Insulin EIA was performed with the following results

Serum Samples

$$y = 0.8061x + 2.8674$$

$$r = 0.997$$

$$n = 32$$

x = another company's Insulin

Plasma Samples

$$y = 0.8161x + 4.5552$$

$$r = 0.986$$

$$n = 47$$

x = another company's Insulin

y = Insulin

Assay Range

1 - 100 μ U/mL

INTERFERENCE:

No cross-reactivity with pro-insulin was observed. Hemoglobin, bile, or rheumatoid factor did not interfere with the assay.

Bilirubin F	No interference up to 19.3 mg/dL
Bilirubin C	No interference up to 19.9 mg/dL
Hemoglobin	No interference up to 450 mg/dL
Lipemia	No interference up to a formazin turbidity of 1,550
RF	No interference up to 450 IU/mL

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CALCIUM ARSENAZO

INTENDED USE

For the quantitative determination of calcium concentration in serum or plasma.

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In plasma, calcium consists of three forms: free, conglutinated with proteins or complex with anions such as phosphate, bicarbonate and citrate. Calcium is an absolutely necessary cation for cell functions.

For example: muscle contraction, bone mineralization, glycogen metabolism, blood coagulation and nerve impulses conduction.

Renal diseases, liver diseases, intestinal malabsorption, acute pancreas inflammation, vitamin D deficiency, adrenal cortical hormone therapy, diuretic treatment and hypoparathyroidism all may result in low levels of total calcium. Hyperparathyroidism, hyperthyroidism, Addison's disease, insusceptible vitamin D or vitamin A excessively, malignant diseases with metastases and sarcoidosis will lead to high levels of total calcium. The ratio is one part sample to 100 parts reagent. The change in absorbance at 650 nanometers is directly proportional to the concentration of calcium in the sample. 2 Arsenazo₄ + Ca²⁺ → Ca²⁺-Arsenazo₄ complex (Blue-purple)

REAGENT CONTENTS

Calcium Reagent contents:

Arsenazo III 2.0 17 mM,

8-Hydroxyquinoline Sulfonate 5.0 mM,

Buffer, Surfactant.

REAGENT PREPARATION

No preparation is required.

REAGENT STORAGE AND STABILITY

Calcium Reagent stored unopened at room temperature stable until the expiration date showed on the bottle label.

DO NOT FREEZE

REAGENT DETERIORATION

The reagent should be discarded if the turbidity has occurred; turbidity may be a sign of contamination.

SPECIMEN COLLECTION AND HANDLING

1. Serum, heparin plasma is suitable for sample. Do not use oxalate, EDTA, or citrate plasma. Whole blood and hemolytic are not recommended for use as a sample. Freshly drawn serum is the preferred specimen.

2. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.

3. Centrifuge samples containing precipitate before performing the assay.

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MANUAL PROCEDURE

1. Label tubes Blank, Standard, Controls, Patients, etc.
2. Transfer 1.0 ml of reagent into each tube.
3. Add 0.01 ml (10 µl) of sample to the respective tubes and mix*.
4. Let stand for at least sixty seconds (60) at room temperature.
5. Zero spectrophotometer with blank at 650 nm
6. Read and record absorbance of all tubes. Final color is stable for sixty minutes (60).

* CHEMISTRY MULTI-PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD

CALCULATIONS

Abs. of Unknown x Conc. of std. = Calcium (mg/dl).

Abs. of Standard

Example: If the absorbance of unknown = 0.74, absorbance of standard = 0.84, concentration of standard = 10 mg/dl, then, $0.74 \times 10 = 8.8$ mg/dl 0.84

NOTE: To convert mg/dl to meq/L, divide mg/dl by two (2).

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UNIT: mg/dl QT: / Q2 / Q3 / Q4: /

PRECISION: 0.1 PC: / ABS: /

Calibration Type: Calibrate + Rgt. Bk Calibration Rule: Two-point linear

INTERFERENCE

1. Hemoglobin levels up to 500 mg/dl, Lipemia levels up to 500mg/dl, Ascorbic acid levels up to 30 mg/dl and Bilirubin levels up to 50 mg/dl were found to exhibit negligible interference.

2. On this method, refer to the work of Young for a review of drug and comprehensive list of substances effect on calcium level.

EXPECTED VALUE

8.4 to 10.2 mg/dl, or 2.1 to 2.6 mmol/L

PRECAUTIONS:

1. For in vitro diagnostic use only.

2. Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

3. Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.

PERFORMANCE CHARACTERISTICS

Linearity: 15.0 mg/dl

Comparison: Studies performed using the present method with a similar method yielded a coefficient of correlation of 0.9950 with a regression equation of $y = 0.9980X - 0.0359$. Sample values from 1.8 to 16.1 mg/dl (N=81).

Precision: Within Run precision for Calcium Reagent Set was determined following a modification of NCCLS EPS-A. Two commercial human serum were assayed for 25 times.

Run-Day precision for Calcium Reagent was determined following a modification of NCCLS EPS-A. Two commercial human serum were assayed five times per day for five days for the total of 25 values.

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INTENDED USE

For the direct, colorimetric determination of calcium in human serum or urine.

INTRODUCTION

More than 99% of body calcium exists in bones and teeth. The remaining 1% is present in blood and soft tissues and serves as a cofactor in blood coagulation, metabolism, and neuromuscular physiology. Serum calcium is present in three different forms: 1) nearly 45% is bound by serum proteins, 2) about 5% is complexed in a non-ionized form and 3) the remaining 50% serum calcium is in an ionic (free) form. It is the physiologically active ionic fraction that is important in terms of biological function.

Many factors influence serum calcium levels: hypercalcemia (increased serum calcium) is observed in hyperparathyroidism, hypervitaminosis, sarcoidosis, myeloma, and certain cancers of the bone. Hypocalcemia (decreased serum calcium) is encountered in hypoparathyroidism, rickets, nephrosis, nephritis, steatorrhea, and pancreatitis. Any decrease in serum proteins frequently results in a decrease of the total serum calcium level. Similarly, an increase in protein such as in myeloma may increase the total serum calcium level. There also appears to be a reciprocal relationship between calcium and phosphorus. Increases in serum inorganic phosphorus are associated with a decrease in serum calcium.

Earlier procedures for the determination of calcium involved precipitation of calcium and subsequent determination of the anion of the precipitating agent. More recently, calcium compounds have been determined by atomic absorption spectrophotometry, which has subsequently been recommended as the reference method for determining total serum calcium. Atomic absorption spectrophotometry involves the use of an expensive and dedicated instrument. With the development of chelating reagents and metallochromic indicators, the atomic absorption methods were rapidly replaced by complex metric procedures, which can measure calcium in the serum directly.

PRINCIPLE

Calcium + O-Cresolphthalein Complexone $\xrightarrow{\text{Alkaline}}$
Medium

Calcium - Cresolphthalein Complexone Complex (purple color)
Calcium reacts with cresolphthalein complexone in 8-hydroxyquinoline to form a colored complex (purple color) that absorbs at 570 nm (550 - 580 nm). The intensity of the color is proportional to the calcium concentration. Color intensifiers and a stabilizer are present to minimize interference by other metallic ions.

REAGENT COMPOSITION

When reconstituted as directed, the reagent for calcium contains the following:

1. **Calcium Color Reagent (A):** O-Cresolphthalein Complexone 0.14 mM, 8-Hydroxyquinoline 13 mM.
2. **Calcium Buffer:** Diethylamide 363 mM, Potassium Cyanide, 2 mM, Non-reactive ingredients, and stabilizers in both reagents A and B.
3. **Calcium Standard:** Calcium Carbonate in dilute hydrochloric acid (10 mg/dl).

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.

CAUTION: *In vitro* diagnostic reagents may be hazardous. Handle in accordance with good laboratory procedures, which dictate avoiding ingestion and eye or skin contact.

2. Reagent (A) and (B) may be irritating to skin. Avoid contact.
3. Reagent (B) contains cyanide and should NOT BE PIPETED BY MOUTH.

REAGENT PREPARATION

1. Combine equal volumes of Calcium Color Reagent (A) and Calcium Buffer (B), mix and let stand for ten (10) minutes at room temperature before use.
2. Reagents should be combined in clean plastic vessels. Water and Glassware containing calcium will react with the reagent. All glassware should be rinsed in diluted hydrochloric acid before use.

REAGENT STORAGE AND STABILITY

1. All reagents should be stored at room temperature (15 - 30°C).

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2. Combined reagent (A and B) is stable for two (2) weeks refrigerated and one (1) week at room temperature. Keep bottles tightly capped to prevent evaporation.

REAGENT DETERIORATION

The reagent should be discarded if:

1. Turbidity has occurred; turbidity may be a sign of contamination.
2. Anti coagulants other than heparin should not be used.
3. The reagent fails to meet linearity claims or fails to recover control values in the stated range.

SPECIMEN COLLECTION

Serum

1. Fasting non-hemolyzed serum is specimen of choice.
2. Anti coagulants other than heparin should not be used.
3. Remove serum from clot as soon as possible since red cells can absorb calcium.
4. Older serum specimens containing visible precipitate should not be used.
5. Tubes with cork stoppers should not be used.
6. Serum calcium is stable for twenty-four (24) hours at room temperature (15 - 30°C), one (1) week refrigerated (2 - 8°C) and up to five (5) months frozen and protected from evaporation.

Urine

1. Collect 24 hours urine in a dry clean container containing 20-30 ml of 6N HCl.
2. Alternatively use 1-2 ml of 6N HCl for random sample.

INTERFERING SUBSTANCES

1. Substances that contain calcium or complex calcium should not come in contact with the test specimen. Examples: EDTA, citrate, oxalate, and fluoride.
2. Specimens from patients receiving bromsulphophthalein (BSP) or EDTA should not be used.
3. For a list of substances affecting the accuracy of calcium values with this procedure refer to the references.

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1. Accurate pipetting devices
2. Test tubes/rack
3. Timer
4. Spectrophotometer able to read at 570 nm

GENERAL INSTRUCTIONS

The reagent for Calcium is intended for use as either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

MANUAL PROCEDURE

1. Prepare working reagent. See "REAGENT PREPARATION".
2. Label tubes Blank, Standard, Controls, Patients, etc.
3. Transfer 1.0 ml of working reagent into each tube.
4. Add 0.02 ml (20 µl) of sample to the respective tubes and mix.
5. Let stand for at least sixty seconds (60) at room temperature.
6. Zero spectrophotometer with blank at 570 nm. (Wavelength range: 550 - 600 nm)
7. Read and record absorbances of all tubes. Final color is stable for twenty minutes (20).

* HTI MULTI-PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

* ALTERNATIVE VOLUMES: (0.05 ml sample to 3.0 ml reagent).

LIMITATIONS

1. The reagent is linear to 20 mg/dl. Samples with values above 20 mg/dl should be diluted 1:1 with saline, re-assayed and the result multiplied by two (2).
2. Lipemic or hemolyzed samples require a serum blank. To prepare a serum blank add 0.05 ml (50 µl) of sample to 3.0 ml distilled water. Mix and read against water at 570 nm. Subtract the absorbance reading from the test reading and perform calculation.
3. Contamination of glassware with calcium (usually from detergents) will adversely affect the test. Use acid-washed glassware or plastic tubes.

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CALCULATIONS

$\text{Abs. of Unknown} \times \text{Conc. of std.} = \text{Calcium (mg/dl)}$
 Abs. of Standard

Example: If the absorbance of unknown = 0.74, absorbance of standard = 0.84, concentration of standard = 10 mg/dl, then,
 $0.74 \times 10 = 8.8 \text{ mg/dl}$

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QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established calcium values may be routinely used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction, or procedural errors.

EXPECTED VALUES

8.5 - 10.5 mg/dl

Children under 12, usually have high normal values, which decrease with aging. It is strongly recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

PERFORMANCE CHARACTERISTICS

1. Linearity: 20 mg/dl

2. Comparison: A study performed with a similar method yield a correlation coefficient of 0.97 with a regression equation of
 $y = 0.94x + 0.53$

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Potassium Reagent Set Калий

Intended Use

For the colorimetric determination of potassium in human serum and plasma.

INTRODUCTION

Potassium is the principle cation of the intracellular fluid. It is also an important constituent of the extracellular fluid due to its influence on muscle activity. Its intracellular function parallels that of its extracellular function, namely influencing acid-base balance and osmotic pressure, including water retention.

Elevated potassium levels (hyperkalemia) are often associated with renal failure, dehydration shock or adrenal insufficiency. Decreased potassium levels (hypokalemia) are associated with malnutrition, negative nitrogen balance, gastrointestinal fluid losses and hyperactivity of the adrenal cortex.^{1,2} In previously described colorimetric methods for determination of potassium or sodium, prior deproteinization of serum or plasma specimen was required. Our improved method is the direct spectrophotometric measurement of potassium in blood or plasma.

PRINCIPLE

The amount of potassium is determined by using sodium tetraphenylboron in a specifically prepared mixture to produce a colloidal suspension.³ The turbidity of which is proportional to potassium concentration in the range of 2 - 7 mEq/L.

REAGENT COMPONENTS

1. Potassium Reagent: Sodium Tetraphenylboron 2.1 mM, preservatives and thickening agents.
2. Potassium Standard: Equivalent to 4 mEq/L.

WARNING AND CAUTION

1. Potassium Reagent Set is for "in vitro diagnostic use" only.
2. Sodium Tetraphenylboron is a corrosive substance. Avoid skin contact or ingestion. DO NOT PIPET BY MOUTH. Flush with water if contact occurs.

STORAGE AND STABILITY

Both reagents are stored at room temperature (18-30°C). The reagents are stable until expiration date indicated on the package label.

REAGENT DETERIORATION

Do not use if:

1. The reagent is very cloudy.
2. The reagent fails to achieve assigned value on fresh control serum.

Specimen collection and storage

1. Serum is recommended.
2. Potassium in serum is stable for at least 2 weeks at 2 - 8°C.
3. Specimens for serum potassium analysis should be free from hemolysis since the high concentration of potassium released from red cells significantly increase the serum levels and this invalidates the test results. Blood specimens should also be separated from the red cells shortly after collection to prevent any leakage of potassium from the intracellular into the extracellular fluid. Plasma from anticoagulants not containing potassium is also suitable.

INTERFERENCES

Turbid or icteric samples produce falsely elevated results. Bilirubin above 40 mg/dl and Urea Nitrogen above 80 mg/dl will produce elevated results. Hemolyzed sera produce elevated results. Sera containing high levels of ammonia should be avoided.

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1. Spectrophotometer.
2. Test tubes/rack.
3. Timer.

PROCEDURE

1. Label test tubes: standard, control, patients, etc. A blank is necessary.
2. Pipette 1.0 mL of Potassium Reagent to all tubes.
3. Add 0.01 mL (10 µl) of samples to respective tubes. Mix and let sit at room temperature for 3 minutes.
4. After 3 minutes, set the wavelength of spectrophotometer to 500 nm, zero spectrophotometer with reagent blank. Read and record the absorbance of all tubes.

* TC - MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

Note: If the spectrophotometer being used requires 2.5 mL reagent, use 0.02 mL (20 µl) of sample to 2.5 mL of reagent. Perform the test as described above.

LIMITATIONS

Our method has been found to be linear between 2-7 mEq/L. It is important to note that our method may not produce accurate results when used with potassium calibrator other than that provided by us. Other products contain preservatives that interfere with this procedure and tend to produce false elevated results. Samples with values above 7 mEq/L should be diluted 1:1 with normal saline, re-assayed and results multiplied by two.

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Abs. of STD

Example: If the absorbance of the unknown = 0.200, the absorbance of the

standard is 0.160 and standard concentration is 4 mEq/L, then

$0.200 \times 4 = 5 \text{ mEq/L}$

0.160

QUALITY CONTROL

Serum controls with known normal and abnormal values should be run

routinely to monitor the validity of the reaction.

EXPECTED VALUES⁷

3.4-5.3 mEq/L

It is strongly recommended that each laboratory establish its own

normal range

PERFORMANCE

1. Linearity: 2-7 mEq/L

2. Sensitivity: Based on an instrument resolution of $A = 0.001$, the present method has a sensitivity of 0.006 mEq/L.

3. Comparison: A comparison study performed between our method

and a similar method resulted in a correlation coefficient of 0.99 with a regression equation of $Y = 1.06X - 0.37$.

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INTENDED USE

Kinetic and endpoint methods for the quantitative determination of creatinine in serum.

INTRODUCTION

Creatinine, an anhydride of creatine, is a waste product formed by the spontaneous dehydration of kidneys.¹ Most of the creatinine is found in muscle tissue where it is present as creatine phosphate and serves as a high energy storage reservoir for conversion to ATP. Independent of diet, serum creatinine concentrations depends almost entirely upon its excretion rate by the kidneys. For this reason, its elevation is highly specific for kidney diseases.² The assay of creatinine has been based on the reaction of creatinine with alkaline picric acid as described by Jaffe. Further modifications have developed the Jaffe reaction into a kinetic assay that is fast, simple and avoids interferences. In the endpoint method, acetic acid is used to destroy the creatinine picric acid complex, resulting in a loss of color, the noncreatinine serum constituents retain their picric acid derived colors and thus the differences in absorbances gives the creatinine concentration.

PRINCIPLE

Alkali
Creatinine + Sodium Picric Acid $\longrightarrow$ Creatinine - Picric acid complex (yellow-orange)
Creatinine reacts with picric acid in alkaline conditions to form a color complex which absorbs at 510 nm. The rate of formation of color is proportional to the creatinine concentration in the sample. In the endpoint method the difference in absorbance measurements after color formation and after acidification yields a creatinine value corrected for interfering substances.

REAGENTS

1. Creatinine Picric Acid Reagent: a solution containing 10 mM picric acid.
2. Creatinine Buffer Reagent: a solution containing 10 mM sodium borate, 240 mM sodium hydroxide and surfactant.
*Important Note: If the Creatinine Buffer Reagent has been subjected to cold temperatures, a white precipitate may form. Warm reagent to 37°C with agitation to dissolve all the precipitate before use.
3. Creatinine standard (5 mg/dl): A solution containing creatinine in hydrochloric acid with preservative.
4. Acetic acid reagent: glacial acetic acid

PRECAUTIONS

1. Creatinine Picric Acid Reagent is a strong oxidizing agent. Avoid contact with skin. WIPE ANY SPILLAGE, SINCE PICRIC ACID IS EXPLOSIVE.
2. Creatinine buffer reagent is an alkali. Acetic acid is acid. Avoid ingestion and contact.

REAGENT PREPARATION

Combine equal volumes of Creatinine Picric Acid Reagent and Creatinine Buffer Reagent, mix well.

REAGENT STORAGE

1. Both reagents are stored at room temperature (18 - 25°C).
2. Combined (working) reagent is stable for up to one (1) week.

REAGENT DETERIORATION

The reagent should be discarded if:

1. Turbidity has occurred, turbidity may be a sign of contamination.
 2. The reagent fails to meet linearity claims or fails to recover control values in the stated range.
- SPECIMEN COLLECTION AND STORAGE**
1. Serum is recommended.
 2. Creatinine in serum is stable for 24 (24) hours at refrigerated temperatures (2 - 8°C) and several months when

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frozen (-20° C) and protected from evaporation and contamination.
3. 24-hour urine specimens must be preserved with 15 grams of boric acid.

INTERFERENCES

A number of substances affect the accuracy of creatinine determination. See Young et al⁵ for a comprehensive list.

MATERIALS PROVIDED

1. Creatinine Picric Acid Reagent.
2. Creatinine Buffer Reagent.
3. Creatinine Standard.
4. Acetic Acid Reagent.

MATERIALS REQUIRED BUT NOT INCLUDED

1. Pipetting devices.
2. Timer.
3. Heating bath.
4. Test tubes/rack.
5. Vessel for combining reagents (glass or plastic).
6. Spectrophotometer with a temperature controlled cuvette.

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KINETIC PROCEDURE

1. Label test tubes: Blank, Standard, Control, Patient, etc.
2. Set the spectrophotometer cuvette temperature to 37°C.
3. Pipette 1.0 ml of working reagent into test tubes.
4. Zero spectrophotometer with the reagent blank at 510 nm. (Wavelength range: 500-520 nm).
5. Add 0.5 ml (50 µl) of sample to reagent mix and immediately place into cuvette.
6. After exactly thirty (30) seconds read and record the absorbance (A₁).
7. At exactly sixty (60) seconds after the A₁ reading, read and record the absorbance (A₂).
8. Calculate the change in absorbance (ΔAbs/min.) by subtracting (A₂ - A₁). See "Calculations".

* HTI - MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD

ALTERNATE VOLUMES

If the spectrophotometer in use requires a volume greater than 1.0 ml for accurate reading, use 0.2 ml (200 µl) sample to 3.0 ml reagent. Perform as above.

CALCULATIONS

The creatinine value of the unknown is determined by comparing its absorbance change with that of a known standard.

$$\Delta Abs. (Unknown) \times \text{Concentration of Standard} = mg/dl$$

$$\Delta Abs. (Standard)$$

Where:

ΔAbs. = Absorbance change between readings (A₂-A₁)

SAMPLE CALCULATION

If ΔAbs. (Unknown) = 0.020
ΔAbs. (Standard) = 0.050
Conc. of Standard = 5 mg/dl

Then:

$$0.020 \times 5 = 2.0 \text{ mg/dl creatinine}$$

KINETIC PROCEDURE LIMITATIONS

Albumin at a concentration of 10.0 gm/dl contributes 0.2 mg/dl to the creatinine value, moderate hemolysis (0.2 gm/dl Hgb), grossly icteric and lipemic samples will give elevated results. Acetoacetate above

10 mg/dl will interfere with the results.
MANUAL ENDPOINT PROCEDURE

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1. Combine equal volumes of Creatinine Picric Acid Reagent and Creatinine Buffer Reagent, mix well.
 2. Label test vial, reagent blank, standard, control and unknown test tubes.
 3. Pipette 3.0 ml of working reagent into test tubes.
 4. Transfer 0.1 ml (100 µl) of sample to its respective tube, distilled water to reagent blank and mix.
 5. Place all tubes in 37°C heating bath for fifteen (15) minutes.
 6. Set wavelength of the spectrophotometer at 510 ± 5 nm and zero the instrument with the reagent blank. Read and record the absorbance of all tubes. These absorbances are the A1 readings.
 7. Add 0.1 ml (100 µl) of acetic acid reagent to each tube, mix and let stand for five (5) minutes at room temperature.
 8. Re-zero the instrument with reagent blank. Read and record absorbances of all tubes. These absorbances are the A2 readings.
- * HTI – MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

CALCULATION OF RESULTS

Use the absorbance readings of the standard and unknown to calculate creatinine values as follows:

$$\frac{(A1 - A2)_{\text{Unknown}}}{(A1 - A2)_{\text{Standard}}} = \text{conc. of standard} = \text{creatinine (mg/dl) in unknown}$$

EXAMPLE OF CALCULATION

Assume the value of the standard to be 5.0 mg/dl and that it gave an A1 of 0.23 and A2 of 0.0 respectively, while unknown gave an A1 of 0.28 and A2 of 0.18. The creatinine concentration of the unknown may then be calculated as follows:

$$\frac{0.28 - 0.18}{0.23 - 0.00} \times 5.0 \text{ mg/dl} = \frac{0.10}{0.23} \times 5.0 \text{ mg/dl} = 2.2 \text{ mg/dl}$$

MANUAL ENDPOINT PROCEDURE NOTE

1. The A2 absorbance reading of the standard should be 0.0 because the addition of acid lowers the pH and destroys all of the creatinine picric complex.
2. The final color development is stable for 30 minutes at controlled room temperature (20 - 25°C).

CALIBRATION

Use the aqueous standard provided.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of normal and abnormal control sera with known creatinine values.

EXPECTED VALUES

Serum: Male 0.9 - 5 mg/dl
Female 0.7 - 1.37 mg/dl

PERFORMANCE CHARACTERISTICS

1. Linearity: 25 mg/dl

2. Comparison: A study performed between this procedure and a similar kinetic procedure yielded a correlation coefficient of 0.99 with a regression equation of $y = 0.96x + 0.06$. Serum and control samples used in the study had creatinine values ranging from 0.9 to 8.3 mg/dl.

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INTENDED USE

For the quantitative determination of Creatine Kinase in human serum.


INTRODUCTION

Creatine Kinase (CK) plays an important role in the energy-storing mechanism of tissue by catalyzing the reversible reaction between creatine and ATP to form creatine phosphate and ADP. CK is distributed in various organs; the highest activities (in decreasing order) are skeletal muscle, heart, and brain.¹ Thus, determination of CK is an aid in diagnosing muscular dystrophy and other diseases of the skeletal muscles, myocardial infarction, hypothyroidism, renal diseases, and/or dysfunction.²

The early procedure for determining CK was based on the rate of ATP formation. A modified method was described by Nielson by adding a sulfhydryl compound and AMP to assure maximum CK activity and inhibit adenylate kinase activity.⁴ Optimized conditions for measuring CK were published by Szasz in 1976 as well as by the Scandinavian committee on enzyme.^{5,6} The above procedure was modified again in 1979 to include EDTA.⁷ The present reagent is a modification of the above revision.

PRINCIPLE

Reaction Sequence:

CK catalyzes the conversion of creatine phosphate and ADP to creatine and ATP. The ATP and glucose are converted to ADP and glucose-6-phosphate by hexokinase (HK). Glucose-6-phosphate dehydrogenase (G-6-PDH) oxidizes the D-glucose-6-phosphate and reduces the nicotinamide adenine dinucleotide (NAD). The rate of NADH formation, measured at 340nm, is directly proportional to serum CK activity.

REAGENT COMPOSITION

When reconstituted as directed, the reagent for CK contains the following:

D-Glucose	20mM
Magnesium++	10mM
Adenosine-5'-Monophosphate (AMP)	50mM
N-Acetylthioester (NAC)	20mM
Creatine Phosphate	30mM
Adenosine-5'-Diphosphate (ADP)	2mM
Oxidized Nicotinamide Adenine Dinucleotide Phosphate	2mM
Glucose-6-Phosphate Dehydrogenase (E.C. 1.1.1.49, G-6-PDH)	3,000UL
Hexokinase (E.C. 2.7.1.1, HK)	3,000UL
EDTA	2mM
Buffer	100mM

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use.
- Exercise the normal precautions required for the handling of all laboratory reagents. Pipetting by mouth is not recommended for any laboratory reagent.

REAGENT PREPARATION

Reconstitute with volume of distilled water specified on each vial, swirl to dissolve.

STORAGE AND STABILITY

The reagent should be stored at 2-8°C prior to reconstitution. The reagent may be used until the expiration date indicated on the package label. After reconstitution, the reagent is stable for twentyfour (24) hours at room temperature or twenty-one (21) days refrigerated (2-8°C).

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If reagent appears damp and clumped, deterioration may have occurred and the product should be discarded.

2. Blank Absorbance

If the reconstituted CK REAGENT without added sample has an absorbance greater than 0.70 at 340 nm versus reagent grade water, the reagent is considered to be unsatisfactory for use and should be discarded.

3. Control Assays

Failure to obtain accurate results in the assay of control materials may indicate reagent deterioration.

4. We cannot guarantee the stability of reagents which have been:

a. transferred from their original containers

b. improperly stored prior to or during use

c. contaminated during use

SPECIMEN COLLECTION

Collect whole blood by non-traumatic venipuncture and allow to clot. Centrifuge and remove serum immediately. Serum is reportedly stable for four (4) hours at room temperature, 8-12 hours at 4°C, and 2-3 days when frozen.

Hemolyzed specimens should not be used because of side reactions that may occur due to adenylate kinase, adenosine triphosphate, and glucose-6-phosphate dehydrogenase liberated from red cells.

INTERFERING SUBSTANCES

Certain drugs and medications may affect the activity of CK, see Young et al. 9

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Sample and reagent pipettes, test vials or cuvettes, timer, test tube rack, thermoregulated flowcell, spectrophotometer, control serum.

GENERAL INSTRUCTIONS

The reagent for CK is intended for use either as automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

PROCEDURE (AUTOMATED)

See appropriate instrument application instructions.

PROCEDURE (MANUAL)

1. Reconstitute reagent according to instructions.

2. Pipette 1.0ml of reagent into appropriate tubes and pre-warm at 37°C for four (4) minutes.

3. Zero spectrophotometer with distilled water at 340nm.

4. Add 0.025 ml (25µl) of sample to the reagent, mix, and incubate at 37°C for two (2) minutes.

5. After two minutes, read and record the absorbance. Return the tubes to 37°C. Repeat readings every minute for the next two minutes.

6. Calculate the average absorbance difference per minute (Δ Abs./min.)

7. The Δ Abs./min. multiplied by the factor 6592 (see Calculations) will yield results in IU/L.

8. Samples with values above 1,200 IU/L should be diluted 1:1 with saline, re-assayed, and the results multiplied by two (2).

NOTE: If the spectrophotometer being used requires a final volume greater than 1.0ml for accurate readings, 3ml of reagent, and 0.1ml (100µl) of sample may be used.

If the spectrophotometer being used is equipped with a temperature-controlled cuvette, the reaction mixture may be left in the cuvette while the absorbance readings are taken.

PROCEDURE LIMITATIONS

1. Some inhibitors of CK activity

a. Excessive Mg^{++} , Cl^- , SO_4^{--}

2-

b. Most heavy earth metals, i.e. Zn^{++} , Cu^{++} , Mn^{++}

c. Iodoacetate, and other sulfhydryl binding agents

d. Excess ADP, citrate, fluoride, L-thyroxine

e. Excess uric acid

2. This procedure measures total CK activity irrespective of its tissue or organ of origin.

3. Lower than expected CK values have been reported in samples having high alkaline phosphatase activity.

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$6.22 \times 0.025 = A \text{ Abs./min.} \times 6592^*$

Where: $\Delta \text{Abs./min.}$ = Average absorbance change per minute
TV = Total reaction volume (1.025)
1000 = Conversion of IU/ml to IU/L
 d = Light path in cm (1.0)
 $\square$ = Millimolar absorptivity of NADH (6.22)
SV = Sample volume in ml (0.025)

Example: If your average absorbance change per minute is 0.015, then $0.015 \times 6592 = 98.9 \text{ IU/L}$.

NOTE: If any of the test parameters are changed, a new factor has to be determined using the above formula.
SI UNITS: To convert to SI Units (nkat/L) multiply IU/L by 16.67.
* If 3ml of reagent and 0.1ml (100ul) of samples are used, then IU/L = $\Delta \text{Abs./min.} \times 4984$.

QUALITY CONTROL

Use control sera with known normal and abnormal values to monitor the integrity of the reaction. Values should be those acceptable for this method and temperature.

TEMPERATURE CONVERSION FACTORS

You may convert results to approximate results at other temperatures by multiplying by the appropriate temperature factor.

Performance Temperature	Reporting Temperature		
	25°C (TF)	30°C (TF)	37°C (TF)
25°C	1.00	1.38	2.30
30°C	0.72	1.00	1.67
37°C	0.4	0.60	1.00

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PERFORMANCE

1. Linearity: 1,200 IU/L

2. Comparison: Studies done between this procedure and a similar procedure yield a correlation coefficient of 0.991 with a regression equation of $Y = 1.01X - 0.29$.

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INTENDED USE


CK-MB is intended to measure the activity of isoenzyme CK-MB in human serum.

INTRODUCTION

Creatine kinases are dimeric molecules composed of M and B subunits and exist as the isoenzymes MM, MB, and BB. The subunits M and B are immunologically distinct; CK-MM and CK-MB are distributed primarily in the skeletal muscle and heart muscle, respectively. While CK-BB is present mainly in the brain and in tissues composed of smooth muscle. 2 Following acute myocardial infarction, CK-MB activity increases significantly and this elevation is highly specific for the laboratory diagnosis of myocardial infarction. 3,4 Although total CK activity usually increases following myocardial infarction, in some patients only the CK-MB activity increases, while the total CK remains in the normal range. 5 In conventional methods, CK isoenzymes are quantitated after first separating the three species by electrophoresis, 6 column anion exchange, 7, or batch anion exchange chromatography. However, these methods are time consuming. In recent times, Wurzburg et al. has introduced an immunoinhibition method. 8 This methodology forms the basis of our CK-MB reagent.

PRINCIPLE

The sample is incubated in the CK-MB reagent which includes the anti-CK-MB antibody. The activity of the noninhibited CK-B is then determined using the following series of reactions:

CK-B catalyzes the reversible phosphorylation of ADP, in the presence of creatine phosphate, to form ATP and creatine. The auxiliary enzyme hexokinase (HK) catalyzes the phosphorylation of glucose by the ATP format, to produce ADP and glucose-6-phosphate (G-6-P) is oxidized to 6-phosphogluconate with the concomitant production of NADH. The rate of NADH formation, measured at 340 nm, is directly proportional to serum CK-B activity.

REAGENT COMPOSITION

1. CK-MB Reagent:
Creatine Phosphate 30 mM, Adenosine-5'-Phosphate 2 mM, Nicotinamide Adenine Dinucleotide (NAD) 2 mM, Hexokinase (Yeast) ≥ 3000 U/L, Glucose-6-Phosphate Dehydrogenase (Bacterial) ≥ 2000 U/L.
2. CK-MB Diluent:
Buffer 100 mM, Anti-Human CK-MB Antibody (Goat)-sufficient amount to inhibit up to 1500 U/L of CK-MM at 37°C.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Exercise the normal precautions required for the handling of all laboratory reagents. Pipetting by mouth is not recommended for any laboratory reagent.

REAGENT PREPARATION

Reconstitute each vial of CK-MB reagent with the volume CK-MB diluent specified on the vial label. Swirl to dissolve.

STORAGE AND STABILITY

The unconstituted reagent and diluent should be stored at 2 - 8°C. They are stable until the expiration date. The reconstituted reagent is stable for at least 7 days in refrigerator (2 - 8°C) and 24 hours at room temperature (15 - 30°C).

SPECIMEN COLLECTION

Serum is the specimen of choice for this assay. Avoid exposure of samples to strong light. Store samples in refrigerator (2-8°C), but no longer than one week. Freezing of samples (-20°C) results in minimal loss of activity.

INTERFERING SUBSTANCES

Extremely hemolyzed samples are not suitable for the test since they may contain high levels of adenylate kinase, ATP, and glucose-6-phosphate, which interfere with the assay to yield false results. Drugs and other substances, which may interfere with the determination of creatine kinase

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activity, have been listed by Young et al 10 The described procedure may overestimate CK-MB values if CK-BB activity in the serum is very high. However, CK-BB activity is usually absent in sera from normal individuals and patients with myocardial infarction¹¹. Some investigators have observed a macro form of BB (immunoglobulin complexed), which may be measured as B in this assay^{12,13,14}. The presence of macro BB in the specimen should be suspected if the CK-B activity measured by this procedure represents more than 20% of the total CK activity.

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Sample and reagent pipettes, test tubes or cuvettes, timer, thermoregulated flowcell, spectrophotometer, control serum.

AUTOMATED PROCEDURE

Consult our appropriate instrument application instructions.

NOTE: Certain instruments require different reconstitution volumes than those stated on the vial label. Refer to appropriate application sheets.

MANUAL PROCEDURE

1. Reconstitute CK-MB reagent according to instructions.
2. Pipette 1.0 ml of CK-MB reagent into the appropriate test tubes and pre-warm at 37°C for at least two (2) minutes.
3. Zero spectrophotometer with water at 340 nm.
4. Add 0.050 ml (50 µl) of sample to the reagent, mix, and incubate at 37°C for five (5) minutes.
5. After five minutes, read and record the change in absorbance per minute for two (2) minutes.
6. Calculate the average absorbance difference per minute ($\Delta Abs./min.$).
7. The $\Delta Abs./min.$ multiplied by the factor 3376 (see Calculations) will yield CK-B results in IU/L.
8. Samples with values above 1500 IU/L should be diluted 1:1 with saline, re-assayed, and the results multiplied by two (2).

NOTE: If the spectrophotometer being used requires a final volume greater than 1.0 ml for accurate readings, 3.0 ml of reagent and 0.15 ml (150 µl) of sample can be used. If the spectrophotometer being used is equipped with a temperature-controlled cuvette, the reaction mixture may be left in the cuvette while readings are taken.

LIMITATIONS

The procedure assumes that CK-BB activity in the sample is negligible. If a significant amount of CK-BB activity is present, then the CK-MB activity will also be overestimated.

CALCULATIONS

1. Total CK Activity:

Determine total CK Activity in serum according to the directions provided in the package insert for CK Reagent.

2. CK-B Activity:

$$IU/L = \frac{\Delta Abs./min. \times TV \times 1000}{d \times \epsilon \times SV} = \frac{Abs./min. \times 1.050 \times 1000}{1 \times 6.22 \times 0.050}$$

$$= \Delta Abs./min. \times 3376$$

Where:

- $\Delta Abs./min.$ = Average absorbance change per minute
- TV = Total reaction volume (1.050)
- 1000 = Conversion of IU/ml to IU/L
- d = Light path in cm (1.0)
- ϵ = Millimolar absorptivity of NADH (6.22)
- SV = Sample volume in ml (0.050)

3. CK-MB Activity:

CK-MB activity is calculated from CK-B activity as follows:

$$CK-MB \text{ Activity (IU/L)} = CK-B \text{ Activity IU/L} \times 2^*$$

*CK-MB molecule is a dimer consisting of a B subunit and an M subunit. Antibody complexing with the M subunit results in loss of half the catalytic activity of the CK-MB molecule. Therefore, CK-MB

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activity in the sample is equal to twice the CK-B activity.

Example:

If your average absorbance change per minute is 0.020, then:
 $0.020 \times 3376 = 67.5$ IU/L (CK-B Activity)

NOTE: CK-MB activity (IU/L) = CK-B activity (IU/L) $\times$ 2

For example, if CK-B activity is 67.5 IU/L, then
CK-MB = $67.5 \times 2 = 135.0$

Percentage of CK-MB activity in sample is:

% CK-MB activity = CK-MB activity $\times$ 100
Total CK activity

For example, if the total CK activity is 1007 IU/L, the CK-B activity is 67.5 IU/L, and the CK-MB activity is 135 IU/L, then

% CK-MB activity = $\frac{135 \times 100}{1007} = 13.5\%$

CALIBRATION

The CK activity in the sample is calculated based on the millimolar absorptivity of NAD. CK-MB reagent is suitable for CK isoenzyme assay when total CK activity in the sample does not exceed 1500 IU/L at 37°C.

QUALITY CONTROL

Use control sera with known normal and abnormal values to monitor the integrity of the reaction in each set of assay. Values should be acceptable for this method and temperature.

TEMPERATURE CONVERSION FACTORS

To convert CK-MB activity at 37°C to 30°C value, multiply the result by 0.60

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PERFORMANCE CHARACTERISTICS

1. **Linearity:** 1500 IU/L

2. **Sensitivity:** Based on an instrument resolution of $A = 0.001$, this procedure has a sensitivity of 4 IU/L.

3. **Correlation:** Studies done between this procedure and Sigma procedure yield a correlation coefficient of 0.98 with a regression equation of $Y = 0.98X - 0.823$ ($N = 40$)

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Intended Use

For the quantitative determination of lactate in human plasma. For *in vitro* diagnostic use only.

Clinical Significance

Lactate determinations are used in the diagnosis of lactate acidosis. Shock is the most widely recognized cause of lactic acidosis although, it is possible for elevated lactate levels to precede shock. Myocardial infarction, severe congestive heart failure, pulmonary edema and blood loss are the common causes of shock which will produce lactic acidosis. Lactic acidosis may also result from renal failure and leukemia. Thiamine deficiency and diabetic ketoacidosis will usually result in increased levels of lactate.

Method History

Originally lactic acid determinations were performed by either titrimetric or colorimetric methods.

The first enzymatic method for lactic acid was based on the transfer of hydrogen from lactate to potassium ferrocyanide by lactate dehydrogenase (LD). This procedure was very cumbersome and did not gain wide acceptance. More current enzymatic methods involved the measurement of NADH formed from the oxidation of lactate by LD 1.2. This method has become more widely used, but still suffers from instability in many analyzer systems. The current enzymatic method is based on the action of lactate oxidase. This method is fast, accurate and is considerably more stable than previous enzymatic methods.

Principle

Lactate oxidase catalyzes the oxidation of lactic acid to pyruvate and hydrogen peroxide. Peroxidase then catalyzes the reaction of hydrogen peroxide with a hydrogen donor, in the presence of 4-aminophenazone, to form a dye. Color intensity, measured at 550nm, is proportional to the lactate concentration in the sample.

Lactate Oxidase

Pyruvate + Hydrogen Peroxide
Peroxidase

Hydrogen Peroxide + TOOS + 4-AAP → Dye (550 nm)

Reagents

Lactate Reagent (R1): TRIS Buffer 100mM, 4-aminopyrene 1.7mM, Peroxidase (Horseradish) > 10,000 U/L, Surfactant, Stabilizer, Sodium Azide (0.09%) as preservative.

Lactate Reagent (R2): TRIS Buffer 100mM, Lactate Oxidase (Microbial) > 1,000 U/L, TOOS 1.5mM, Surfactant, Stabilizer, Sodium Azide (0.09%) as preservative.

Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Reagents contain sodium azide as preservative. Upon disposal flush with large volumes of water.
3. All specimens used in this test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.
4. Do not use the reagents beyond the expiration date printed on the kit label.

Reagent Preparation

Lactate reagents R1 and R2 are ready to use for instruments suitable for two reagent analysis.

Reagent Storage

All reagents are stable until the expiration date on the label when stored at 2-8°C.

Specimen Collection and Storage

Plasma collected in sodium fluoride/potassium oxalate is the recommended specimen. The specimen should be immediately placed on ice and the cells must be separated within 15 minutes. 3 The sample should be drawn from a status-free vein. 4 If not analyzed promptly, specimens may be stored at 2-8°C for up to 2 days. If specimens need to be stored for more than 2 days, they may be stored for one month frozen at -20°C.

Interferences

All interference studies were conducted based on the procedures recommended in NCCLS guideline No. EP7-P.6 Hemoglobin at levels up to 500 mg/dl and Bilirubin at levels up to 20 mg/dl were found to exhibit negligible interference (<5%) on this method. Samples with levels of interfering substances higher than the upper limits should be diluted with physiological saline before assaying.

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Multiply the result obtained from the manual dilution by the appropriate dilution

factor. For a comprehensive review of drug interference on lactate levels see Young et al.7

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1. Lactate standard or suitable serum-based calibrator
2. Controls with normal and elevated levels of lactate.
3. Automated clinical chemistry analyzer capable of accommodating two reagent assays.

Procedure

Below is a general example of the lactate test procedure for an automated analyzer. For assistance with applications on automated analyzers, please contact the manufacturer's Technical Service Department.

Sample + Reagent 1 ————— Reagent 2 ————— Measurement of
Bul 300ul 30 sec. 200ul 5 min absorbance (700/546nm)

Lactate results

Limitations

1. Anticoagulants containing citrate should not be used.
2. Protect the reagents from direct sunlight.
3. Samples with values greater than 20mmol/L must be diluted 1:1 with saline and re-assayed. Multiply the result by two.

Calibration

Use an NIST-traceable lactate standard, or a suitable serum-based lactate standard. The calibration should be calibrated according to the instrument manufacturer's instructions. If control results are found to be out of range, the procedure should be re-calibrated.

Quality Control

Reliability of test results should be routinely monitored with control materials that reasonably emulate performance of patient specimens. Quality control materials are intended for use only as monitors of accuracy and precision. The recovery of control values within the appropriate range should be the criteria used in evaluation of future assay performance. Controls should be run with every working shift in which lactate assays are performed. It is recommended that each laboratory establish their own frequency of control determination. Quality control requirements should be determined in conformance with local, state, and/or Federal regulations or accreditation requirements.

Results

To convert from S.I. units to conventional units, multiply the S.I. units by 9.01.

Example: mmol/L x 9.01 = mg/dL Lactate

Expected Values

The following reference range is suggested for L-Lactate.8

Venous 0.5-2.2 mmol/L

Arterial 0.5-1.6 mmol/L

It is highly recommended that each laboratory establish its own range of expected values.

Performance

1. Assay Range: 0-20 mmol/L.
2. Comparison: This lactate reagent was compared to the method performed on the Dade Chemistry Analyzer. The study was performed using 57 patient samples ranging from 0.3-10.4 mmol/L. Data was subjected to least-squares linear regression analysis which yielded a correlation coefficient of 0.998 with a regression equation of $y = 0.97x + 0.1$.
3. Precision:

Within-Day precision for the Lactate Reagent was

determined following a modification of NCCLS document EP5-T2.9

Within-Day precision studies produced the following results:

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Day To Day precision was also determined following a modification of NCCLS document EPS-T2.9 Day to Day precision studies produced the following results:

Day To Day

Sample	N	Mean	S.D.	C.V.%
Low	20	1.51	0.04	2.65
Mid	20	4.12	0.09	2.18
High	20	9.19	0.17	1.85

4. Sensitivity: The analytical sensitivity for lactate was determined to be 0.15 absorbance units per 1 mmol/L of lactate.

References

1. Gutmann, I., Wahlefeld, A., Methods of Enzymatic Analysis, 2nd Ed., Academic Press, New York, 1974, 1464.
2. Noll, F., Methods of Enzymatic Analysis, 2nd Ed., Academic Press, New York, 1974, 1465.
3. Tietz, N.W., Fundamentals of Clinical Chemistry, 4th Ed., W.B. Saunders Company, Philadelphia, 1996, 367.
4. Tietz, N.W., Clinical Guide to Laboratory Tests, 3rd Ed., W.B. Saunders Company, Philadelphia, 1995, 382-383.
5. Westgard, J.O., Lahmeyer, B.L., Birnbaum, M.L., Clin Chem 1972, 18:1334-1338.
6. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 4, No. 8, June 1984.
7. Young, D.S., effects of Drugs on Clinical Laboratory Tests, 3rd Ed., AACCPress, Washington D.C., 1990.
8. Tietz, N.W., Fundamentals of Clinical Chemistry, 4th Ed., W.B. Saunders Company, 1996, 801.
9. NCCLS document "Evaluation of Precision Performance of Clinical

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HTI-L533-600

INTENDED USE

Lactate Dehydrogenase LS (LDH-L) liquid reagent is intended to use for the *In Vitro* quantitative kinetic determination of Lactate Dehydrogenase activity in serum.

INTRODUCTION

The enzyme lactate dehydrogenase (LDH-L) is distributed in tissues particularly heart, liver, muscle, and kidney. The enzyme found in circulation is a mixture of five isoenzymes based on their mobility. Elevated serum levels of LDH-L are found in serum in myocardial infarction, liver disease, renal disease, certain forms of anemia, malignant diseases, and progressive muscle dystrophy.^{1,2}

LDH catalyzes the following reaction:

The LDH-L enzyme activity can be measured in both directions. Using optimal conditions for both directions and taking into account isoenzyme variations, assays in either direction are considered to be equivalent.³ However, Lactate-to-Pyruvate method offers a number of advantages:⁴ (1) the rate of reaction is linear over a wide range (2) no pre-incubation is required and (3) better reagent stability.

PRINCIPLE

LDH catalyzes the oxidation of lactate to pyruvate in the presence of NAD, which is subsequently reduced to NADH. The rate of NADH formation measured at 340 nm is directly proportional to serum LDH-L activity.

REAGENT COMPOSITION

After combining R1 and R2 as directed, the reagent contains: Lithium lactate 55 mmol/L, NAD 7.5 mmol/L, Buffer (8.95), and stabilizers and preservatives.

PRECAUTIONS

1. For "in vitro" diagnostic use only.
2. The reagent contains sodium azide as a preservative. Do Not Ingest. Avoid skin and eye contact.
3. All specimens and controls should be handled as potentially infectious, use safe laboratory procedures. (NCCLS M29-T2)

REAGENT PREPARATION

The working reagent is prepared by mixing five (5) volumes of R1 with one (1) volume of R2 in a disposable container.

Example: 25ml R1 + 5ml R2

REAGENT STORAGE

1. Store reagent at 2 - 8°C.
2. The reagents are stable until the expiration date if stored as directed.
3. The "working" reagent is stable for 14 days at 2-8°C.
4. Protect from direct light.
5. Avoid microbial contamination.

REAGENT DETERIORATION

1. If the reagent blank before serum addition has an absorbance that exceeds 0.600 at 340 nm, the reagent may have deteriorated.
2. Failure to obtain accurate results in the assay of control materials may indicate reagent deterioration.

SPECIMEN COLLECTION AND HANDLING

1. Serum with any visible hemolysis cannot be used because of the contamination of this sample with large amount of LDH released from the erythrocytes.
Serum should be separated from the clot promptly.
Samples should be assayed soon after collection. LDH in serum is reported stable for two to three days at room temperature.²
2. The liver LDH is particularly labile and is destroyed if frozen and thawed.⁵

INTERFERENCE

1. Oxalate, oxamates, and EDTA will inhibit LDH.
2. Young, et. al. gave a list of drugs and other substances interfere with the determination of LDH activity.⁶

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MATERIALS PROVIDED

Lactate Dehydrogenase LS* Buffer Reagent
Lactate Dehydrogenase LS* Co-Enzyme Reagent

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- 1 Accurate pipetting devices.
- 2 Test tubes and rack.
- 3 Timer.
- 4 Heating bath or block (37°C).
- 5 Spectrophotometer capable of reading at 340 nm (UV).

AUTOMATED PROCEDURE

Refer to appropriate instrument application instruction.

MANUAL PROCEDURE

1. Prepare the "working" reagent as directed.
2. Pipette 1.0 ml of reagent into appropriate tubes and prewarm at 37°C for at least five (5) minutes.
3. Zero spectrophotometer with water at 340 nm.
4. Add 0.025 ml (25 µl) of sample to reagent, mix and incubate at 37°C for one (1) minute.
5. After one (1) minute, read and record the absorbance and return tube to 37°C. Repeat reading every minute for the next two (2) minutes.
6. Calculate the average absorbance difference per minute (ΔAbs/min.)
7. The absorbance ΔAbs/min. multiplied by 6592 will yield results in U/L.

LIMITATIONS

Samples with values above 1,000 U/L should be diluted 1:1 with saline, re-assayed and the results multiplied by two (2) to compensate for the dilution.

CALIBRATION

The procedure is standardized by means of the millimolar absorptivity of NADH, which is 6.22 at 340 nm under the test conditions described.

CALCULATION

Lactate Dehydrogenase activity is expressed as U/L.

Unit definition: One international unit (U/L) is the amount of enzyme that will reduce one micromole of NAD per minute at specific temperature.

$$U/L = \frac{\Delta Abs/min \times TV \times 1000}{d \times e \times SV} = \frac{\Delta Abs/min \times 1.025 \times 1000}{1.0 \times 6.22 \times 0.025}$$

$$= \Delta Abs/min \times 6592$$

Where:

A Abs/min	= Average absorbance change per minute
TV	= Total reaction volume (1.025)
1000	= Conversion of U/ml to U/L
d	= Light path in cm (1.0 cm)
e	= Millimolar absorptivity of NADH (6.22)
SV	= Sample volume in ml (0.025)

Example:

If the average absorbance change per minute is 0.022, then $0.022 \times 6592 = 145$ U/L.

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with LDH-L values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction or procedural errors.

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PERFORMANCE CHARACTERISTICS

1. Sensitivity: Based on an instrument resolution of 0.001 absorbance, this procedure has a sensitivity of 3.4 IU/L.

2. Linearity: Linear from 0 - 1,000 IU/L.

3. Comparison: A group of 54 sera ranging in Lactate Dehydrogenase values from 15 - 992 IU/L was assayed by this method and by a similar commercially available reagent. Comparison of the results yielded a correlation coefficient of 0.996 and the regression equation was $y = 0.954x - 9.02$.

(Comparison studies were Tentative Guideline, EP9-T.)

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REFERENCES

- Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods W.B. Saunders and Company, Philadelphia, PA p. 365 (1979).
Tietz, R.W., Fundamentals of Clinical Chemistry, W.B. Saunders and Company, Philadelphia, PA p. 652 (1976).
Henry, R.J. et al., Clinical Chemistry Principles and Techniques, 2nd Ed., Harper and Row, Hagerstown (MD) p. 819-831 (1974).
Amador, E., et al. Clin. Chem. 9:391 (1963).
Kreutzer, H.H., et al. Clin. Chem. Acta 9:64 (1964).
Young, D. S., et al. Clin. Chem., 21: ID 432D (1975).

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INTENDED USE

Lipase reagent is used for the quantitative turbidimetric determination of pancreatic lipase in human serum.

INTRODUCTION

Lipase is defined as that group of enzymes, which hydrolyze the glycerol esters of long-chain fatty acids. The measurement of lipase activity in serum and other fluid is to evaluate conditions associated with pancreas. Voget et al. proposed an olive oil emulsion in measuring the rate of change in turbidity over a specific unit of time. Later, Shihabi et al. modified the previous method and eliminated some of interference. Our method is based on the above modifications.

PRINCIPLE

Triglyceride + H₂O → mono- + di-glycerides + fatty acids

Serum lipase hydrolyzes the olive oil emulsion. The decrease in turbidity at 400 nm (after incubation) is proportional to lipase activity in the specimen.

REAGENT (MATERIALS PROVIDED)

1. Substrate: 0.8% (w/v) olive oil in alcohol
2. Buffer: Tris Buffer 69 mM, Sodium Deoxycholate 10 mM, pH 9.0 (37°C).

WARNINGS AND PRECAUTIONS

1. This reagent is for "in vitro" diagnostic use only.
2. Exercise the normal precautions required for the handling of all laboratory reagents. Pipetting by mouth is not recommended for any laboratory reagent.

REAGENT PREPARATION

1. Add lipase buffer to a 50 ml Erlenmeyer flask. Add 25 ml distilled water and swirl to dissolve.
2. Pipette 1 ml of well-mixed lipase substrate into buffer solution. Note: The absorbance of the emulsion prior to use must be greater than 1.0. Due to variations in regional temperatures, the absorbance may be less than 1.0. If this occurs, add 0.5 – 1.0 ml more substrate until absorbance is greater than 1.0.

REAGENT STORAGE

1. Unreconstituted reagent should be stored at room temperature (25 - 30°C).
2. Reconstituted reagent is stable for seven days refrigerated (2 - 8°C) and tightly capped.

REAGENT DETERIORATION

Do not use the emulsion if the absorbance is more than 0.3 below the absorbance of the emulsion when fresh. Prepare a fresh emulsion.

SPECIMEN COLLECTION STORAGE

1. Use fresh serum specimens.
2. Lipase activity in serum is stable at room temperature for one week; sera may be stored for three weeks in the refrigerator (4-8°C) and for several months if frozen. Bacterial contamination of the specimens may result in an increase in lipase activity.

INTERFERING SUBSTANCES

1. Hemolyzed specimens should not be used.
2. A number of drugs and substances affect lipase activity. For a comprehensive list, see Young, et al.

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1. Accurate pipetting devices (3.0 ml and 0.1 ml).
2. 25 ml measuring vessel.
3. 50 ml Erlenmeyer flask.
4. Timer.
5. Test tubes/rack.
6. Spectrophotometer with temperature controlled cuvette.
7. Heating bath (37°C).
8. Serum controls.

PROCEDURE (MANUAL)

1. Reconstitute lipase reagent according to instructions.
2. Label test tubes "blank", "control", "patient", etc.
3. Pipette 3.0 ml of reagent into appropriate test tubes and pre-warm at 37°C for at least five (5) minutes.
4. Zero spectrophotometer with distilled water at 400 nm (Wavelength range: 390-420 nm).

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5. Read and record absorbance of blank and place back in heating bath.
6. Using timed intervals, add 0.1 ml (100 µl) of sample to each tube, mix, and read initial absorbance. Return each tube to heating bath after initial reading.
7. Exactly five (5) minutes after the initial absorbance reading, using the same timed intervals, remove each tube from the heating bath and mix each tube. Read the absorbance of the blank and each sample tube against distilled water.

PROCEDURE NOTES

1. If the Abs. of the "blank" is a negative value, consider it zero.
2. Elevated blank rates (i.e. (0.005) and above may be caused by olive oil coating on cuvette surface. Periodically rinse with acetone followed by water flush.
3. Turbid samples should be diluted with distilled water (1:5). Multiply final answer by dilution factor.
4. Use fresh sera, when possible, for greatest accuracy.

CALIBRATION

The lipase activity in the sample is calculated based on the millimolar absorptivity of olive oil (3.15 in working solution).

QUALITY CONTROL

1. Fresh control sera with known normal and abnormal values should be run for each day to monitor the validity of the reaction.
2. Serum lipase activity should be valid by titrimetric method.

CALCULATIONS

The enzyme activity is expressed in International Units. One International Unit (IU) is defined as the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under decreed conditions.

$$\text{Units/Liter} = \frac{\text{Corrected } \Delta \text{Abs./min}}{\text{Initial blank absorbance}} \times 1953$$

$$\text{Corrected } \Delta \text{Abs./min} = \Delta \text{Abs test} - \Delta \text{Abs blank}$$

Example:

(blank) Initial Abs.	= 0.970
(blank) 5 min. Abs	= 0.967 (if negative result, treat as zero)
$\Delta \text{Abs blank}$	= $0.970 - 0.967 = 0.003$
(test) Initial Abs.	= 1.300
(test) 5 min. Abs	= 1.271
$\Delta \text{Abs test}$	= $1.300 - 1.271 = 0.029$
Corrected $\Delta \text{Abs./min}$	= $0.029 - 0.003 = 0.026$
Therefore:	$\frac{0.026}{0.026} \times 1953 = 52 \text{ IU/L}$

DERIVATION OF FACTOR (1953)

$$\frac{\text{Corrected } \Delta \text{Abs./min} \times 315 \times 3.1}{\text{Initial Abs/blank} \times 0.1} = \frac{315 \times 3.1}{0.1} = \frac{9765}{\Delta \text{Abs./min} \times 5} = 1953$$

315 = concentration of olive oil (micromoles/liter) in the working solution.

3.1 = volume of reaction mixture.

0.1 = volume of sample in ml

($\Delta \text{Abs./min}$) $\times 5$ = conversion to $\Delta \text{Abs./min}$

Note:

To convert IU/L to Cherry-Crandall units, divide IU/L by 70.

Example:

$$\frac{52 \text{ IU/L}}{70} = 0.74 \text{ Cherry-Crandall units}$$

LIMITATIONS

1. Samples with values above 280 IU/L should be diluted 1:1 with distilled water, re-assayed and final results multiplied by two.
2. Bacterial contamination of the specimens may result in an increase in lipase activity.

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INTENDED USE

For the direct quantitative determination of low density lipoprotein cholesterol (LDL-C) in human serum or plasma on automated analyzer For *in vitro* diagnostic use only.

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Plasma lipoproteins are spherical particles that contain varying amounts of cholesterol, triglycerides, phospholipids, and proteins. The phospholipid, free cholesterol and protein constitute the outer surface of the lipoprotein particle, the inner core contains mostly esterified cholesterol and triglycerides. These particles serve to solubilize and transport cholesterol and triglycerides in the bloodstream.

The relative proportions of protein and lipid determine the density of these plasma lipoproteins and provide a basis for their classification. The classes are: very low density lipoproteins (VLDL), low density lipoproteins (LDL), and high density lipoprotein (HDL). Numerous clinical studies have shown that the different lipoprotein classes have varied effects.²⁴ The studies all point to LDL cholesterol as the key factor in the pathogenesis of atherosclerosis and coronary artery disease (CAD), while HDL cholesterol has often been observed to have a protective effect. Even within the normal range of total cholesterol concentrations, an increase in LDL cholesterol can occur with an associated risk for CAD.²⁴

Over the years a variety of methods have been employed for the determination, or estimation, of LDL cholesterol. The Friedewald equation, in a variety of forms, has been most frequently used for the estimation of LDL cholesterol. However, its usefulness is limited and its accuracy has been questioned. Determination of LDL cholesterol by beta-quantification is recognized as the reference method, but the procedure is so cumbersome relatively few laboratories use this method. A recent method using immunoseparation has become popular. However, this method is still requires sample pre-treatment prior to cholesterol determination, making it unsuitable for full automation of the procedure. The method presented here offers direct determination of LDL cholesterol in a two part, liquid stable reagent that is easily adapted to most automated chemistry analyzers.

Principle

The Direct LDL Cholesterol Reagent is a two-part, liquid stable method for directly measuring LDL-C levels in serum or plasma. The method depends on the properties of a unique detergent which eliminates the need for any off-line pre-treatment or centrifugation steps. This detergent (Reagent 1) solubilizes only the non-LDL lipoprotein particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non-color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL-C in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

REAGENT COMPOSITION

1. Direct LDL Cholesterol Reagent 1:

Buffer 100 mmol/L (pH 7.0), Cholesterol esterase from *Pseudomonas*, 800 U/L, Cholesterol oxidase from *Nocardia* sp, 500 U/L, Peroxidase from *Horseradish* 800 U/L, 4-aminoantipyrine, 1 mmol/L, Preservative.

2. Direct LDL Cholesterol Reagent 2:

Buffer 100 mmol/L (pH 7.0), 2'-N,N-bis(4-sulfobutyl)-m-Toluidine-disodium (DSBMT) 1.2%, Preservative.

WARNINGS AND PRECAUTIONS

1. Reagent is intended for *in vitro* diagnostic use only.
2. Do not pipette by mouth.
3. All specimens used in this test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.
4. Do not use the reagents beyond the expiration date printed on the kit label.

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STORAGE AND STABILITY

All reagents are stable until the expiration date on the label when stored at 2 to 8°C.

SPECIMEN COLLECTION AND STORAGE

Serum, EDTA-treated or heparinized plasma are the recommended specimens. Patients are not required to fast prior to blood collection.
Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).
Plasma: Specimens may be collected in EDTA or heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).
If not analyzed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens must be stored for more than 5 days, they may be frozen at -80°C.

INTERFERENCES

All interference studies were conducted according to the procedures recommended in NCCLS guideline No. EP7-P for interference testing in clinical chemistry.¹⁰ Hemoglobin at levels up to 400 mg/dl, Bilirubin at levels up to 20 mg/dl and Triglycerides to 1500 mg/dl were found to exhibit negligible interference (<5%) on this method. Samples with levels of interfering substances higher than the upper limits should be diluted with physiological saline before assaying. Multiply the result obtained from the manual dilution by the appropriate dilution factor. For a comprehensive review of drug interference on serum LDL cholesterol levels see Young et al.¹¹

MATERIALS PROVIDED

1. Direct LDL Reagent 1 (ready to use).
2. Direct LDL Reagent 2 (ready to use).

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1. Direct HDL/LDL Calibrator (Cat. H514-A).
2. LDL cholesterol controls.
3. Automated clinical chemistry analyzer capable of accommodating two-reagent assays.

PROCEDURE

Below is general example of the Direct LDL test procedure for an automated analyzer. All analyzer applications should be validated in accordance with NCEP and CLIA recommendations.¹⁰

1. Use 3ul sample with 300ul of Direct LDL Cholesterol Reagent 1.
2. Equilibrate to 37°C for 5 minutes.
3. Add 100ul of Direct LDL Cholesterol Reagent 2.
4. Equilibrate to 37°C for 5 minutes.
5. Measurement (Absorb. Difference between 660nm & 546nm)

LIMITATIONS

1. Anticoagulants containing citrate should not be used.
2. Protect the reagents from direct sunlight.
3. Samples with values greater than 500 mg/dl must be diluted 1:1 with saline and re-assayed. Multiply the result by two.

CALIBRATION

The Direct HDL/LDL Cholesterol Calibrator is required for calibration. Calibrate with each bottle change or lot change or if control results are found to be out of range. The values of the calibrator were assigned by procedures traceable to the National Reference System for Cholesterol (NRS/CHOL). Refer to Direct HDL/LDL Cholesterol Calibrator package insert for instructions. If control results are found to be out of range, the procedure should be re-calibrated.

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QUALITY CONTROL

Reliability of test results should be routinely monitored with control materials that

reasonably emulate the performance of patient specimens.¹⁰ Quality control materials are intended for use only as monitors of accuracy and precision. Any HDL/LDL controls would be suitable for use with this assay. The recovery of control values within the appropriate range should be the criteria used in evaluation of future assay performance. Controls should be run with every working shift in which LDL-C assays are performed. It is recommended that each laboratory establish its own frequency of control determination. Quality control requirements should be determined in conformance with local, state, and/or Federal regulations or accreditation requirements.

RESULTS

Concentration of patient (mg/dl) =

Abs (patient) X Concentration of standard (mg/dl)
Abs (Standard)

Example: Abs (patient)=0.32 Abs (standard)=0.28
Concentration of standard=114 mg/dl

Concentration of patient (mg/dl)= $\frac{0.32}{0.28} \times 114 = 130.2$ mg/dl

To convert from conventional units to S.I. units, multiply the conventional units by 0.02586.

Example: mg/dl x 0.02586 = mmol/L LDL-C

EXPECTED VALUES

The following NCEP recommendations for patient classifications are suggested for the prevention and management of coronary heart disease:

LDL Cholesterol	Classifications
<100mg/dl (2.586mmol/L)	Optimal
100-129mg/dl (2.586-3.34mmol/L)	Near Optimal/Above Optimal
130-159 mg/dl (3.36-4.11mmol/L)	Borderline High Risk
160-189mg/dl (4.14-4.89mmol/L)	High Risk
≥190 mg/dl (4.91mmol/L)	Very High Risk

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PERFORMANCE CHARACTERISTICS

Assay Range: 5-500 mg/dl

Accuracy: Accuracy of the Direct LDL Cholesterol Reagent method was verified by comparison to the reference method (AutoLDL Cholesterol) on Hitachi 717. There were 52 serum samples in this study. The results of this study are below.

Method	Direct LDL(Y) Cholesterol	AutoLDL(X) Cholesterol
N	52	52
Range (mg/dl)	34.5-270.3	30.5-274.4
Regression Analysis	Y=0.94x+3.39 mg/dl	
Correlation Coefficient	R=0.987	

Precision:

Within-Day precision for the Direct LDL Cholesterol Reagent was determined on Hitachi 717 following a modification of NCCLS document EP5-A. Within-Day precision studies produced the following results:

Sample	Sample 1	Sample 2
N	25	25
Mean LDL Cholesterol (mg/dl)	178	33

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Day-to-Day precision was also determined on Hitachi 717 following a modification of NCCLS document EP5-A.¹⁷ Day-to-Day precision studies produced the following results:

	Sample 1	Sample 2
N	25	25
Mean LDL Cholesterol (mg/dl)	174.2	32.7
Standard Deviation (mg/dl)	5.04	0.9
Coefficient of Variation (%)	2.8	2.7

Sensitivity. The analytical sensitivity for Direct LDL Cholesterol was determined by Hitachi 717 as 1 mg/dl of LDL cholesterol.

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REFERENCES

- Goto, A.M., Lipoprotein Metabolism and the etiology of Hyperlipidemia, Hospital practice, 23: Suppl. 1, 4 (1988).
- Crouse, J.R., et al., Studies of Low Density Lipoprotein Molecular Weight in Human Beings with Coronary Artery Disease, J. Lipid Res., 25: 586 (1985).
- Badimon, J.J., Badimon L., Fuster V., Regression of Atherosclerotic Lesions by High-density lipoprotein Plasma fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85: 1234-41 (1990).
- Castelli, W.P., et al., Cholesterol and other Lipids in coronary heart disease, Circulation, 55: 767 (1977).
- Barr, D.P., Russ, E.M., Elder, H.A., Protein-Lipid Relationships in Human Plasma, Am. J. Med. 11: 480 (1951).
- Gordon, T., et al., High Density Lipoprotein as a Protective Factor Against Coronary Heart Disease, Am. J. Med., 62: 707 (1977).
- William, P., Robinson, D., Baily A., High Density Lipoprotein and Coronary Risk Factor, Lancet, 1: 72 (1979).
- Kannel, W.B., Castelli W.P., Gordon, T., Cholesterol in the Prediction of Artherosclerotic Disease, New Perspectives Based on the Framingham Study, Am. Intern. Med., 90: 85 (1979).
- National Institutes on Health Publication no. 93-3095, September 1993.
- Warnick, G. Russell, Wood Peter D., National Cholesterol Education Program Recommendations for Measurement of High Density Lipoprotein Cholesterol, Executive Summary, Clinical Chemistry, Vol. 41, No. 10, 1995.
- Grundy, S.M., et al., Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II), JAMA 1993, 269: 23, 3015-3023.
- National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 4, No. 8, June 1984.
- Young, D.S. Effects of Drugs on clinical Laboratory Tests, 3rd ed., AACCC Press, Washington, D.C., 1990, 3-104 thru 3-106.
- Tietz, N.W., Clinical Guide to Laboratory Tests, W.B. Saunders Co., Philadelphia, 1986, p. 256.
- Carey, R., Gerber, C.C., Evaluation of Methods, In Kaplan LA, Pesce, A.J., eds. Clinical Chemistry: theory, analysis and correlation. Third Edition. St. Louis: The CV Mosby Company.
- Westgard, J.O., Carey, R.N., Wold, S., Criteria for judging precision and accuracy in method development and evaluation. Clinical Chemistry 1974; 20: 825-833.
- NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices" 2nd Ed. 1992.

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htimed.com

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INTENDED USE

For the quantitative determination of high-density lipoprotein (HDL) in human serum or plasma on automated analyzer. For *in vitro* diagnostic use only.

SUMMARY AND EXPLANATION OF THE TEST

Cholesterol is a fatty substance found in blood, bile and brain tissue. It serves as a precursor to bile acids, steroids and vitamin D. The concentration of total cholesterol in serum has been associated with metabolic, infectious and coronary heart diseases. In the plasma, cholesterol is transported by three lipoproteins: high-density lipoprotein (HDL-Cholesterol), low density lipoprotein (LDL-Cholesterol), and very low density lipoprotein (VLDL-Cholesterol).

The role of HDL particles in lipid metabolism is primarily the uptake and transport of cholesterol from peripheral tissue to the liver. This process is known as reverse cholesterol transport and has been proposed as a cardio protective mechanism. Low HDL-C levels have repeatedly been associated with an increased risk of coronary heart disease and coronary artery disease. Thus the determination of serum HDL cholesterol has been recognized as a useful tool in identifying high-risk patients.

The CDC reference method for HDL cholesterol uses ultracentrifugation followed by chemical precipitation to separate HDL from other lipoproteins, followed by cholesterol measurement using a modified Abel-Kendall assay. This method is considered too time consuming and labor intensive for use in routine analysis. Historically, most laboratories have used one of several methods for the selective precipitation and removal of LDL and VLDL, followed by the enzymatic measurement of HDL-C in the supernatant fraction. Since almost all of these methods required manual separation steps, HDL cholesterol determinations could not be fully automated. Also, the dilution of the sample resulted in an enzymatic determination of cholesterol with low sensitivity. As a result, the routine determination of HDL cholesterol has suffered from both long turnaround times and poor reproducibility.

PRINCIPLE

The Direct HDL Cholesterol assay is a homogeneous method for directly measuring serum HDL-C levels without the need for any off-line pretreatment or centrifugation steps. The method is in a two-reagent format. The first reagent stabilizes LDL, VLDL, and chylomicrons. The second reagent contains PEG modified enzymes that selectively react with the cholesterol present in the HDL particles. Consequently, only the HDL cholesterol is subject to cholesterol measurement.

REAGENT COMPOSITION

1. Direct HDL Cholesterol Reagent 1:
Magnesium chloride 100mM, Aminoantipyrene 1 mMol/L, buffer, pH 7.0±0.1, preservative.
2. Direct HDL Cholesterol Reagent 2:
Peroxidase from *Horseradish* (POD) 4KU/L, Cholesterol Oxidase from *Nocardia* sp. (PEG-CO) 1 KU/L, Cholesterol Esterase from *Pseudomonas* (PEG-CE) 1 KU/L, N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (HDAOS) 0.3g/L, buffer, pH 7.0±0.1, surfactant, preservative.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Exercise the normal precautions required for the handling of all laboratory reagents. Pipetting by mouth is not recommended for any laboratory reagent.
3. All specimens used in this test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.
4. Do not use the reagent after the expiration date printed on the kit.

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Store the reagent set at 2 - 8 °C. The reagent is stable until the expiration date indicated on the bottle label.

SPECIMEN COLLECTION AND STORAGE

Serum, EDTA-treated or heparinized plasma are the recommended specimens.
Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).
Plasma: Specimens may be collected in EDTA or heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).
If not analyzed promptly, specimens may be stored at 2-8°C for up to 1 week. If specimens need to be stored for more than 1 week, they may be preserved at less than -20°C for up to 1 month. For storage periods of 1 month to 2 years, samples should be preserved at -70°C.

INTERFERENCES

All interference studies were conducted according to the procedures recommended in NCCLS guideline NO. EP7-P for interference testing in clinical chemistry.¹² Hemoglobin levels up to 100 mg/dl, Triglyceride levels up to 1800 mg/dl and Bilirubin levels up to 20 mg/dl were found to exhibit negligible interference (<5%) on this method. Samples with levels of interfering substances higher than the upper limits should be diluted with physiological saline before assaying. Refer to the work of Young for a review of drug effects on serum HDL cholesterol levels.

MATERIALS PROVIDED

1. Direct HDL Cholesterol Reagent 1 (Ready to use)
2. Direct HDL Cholesterol Reagent 2 (Ready to use)

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1. Direct HDL/LDL Cholesterol Calibrator (Cat. # H514-A)
2. Lipid Control Set
3. Automated clinical chemistry analyzer capable of accommodating two-reagent assays.

PROCEDURE

Below is a general example of the Direct HDL Cholesterol test procedure for an automated analyzer. All analyzer applications should be validated in accordance with NCEP and CLIA recommendations.¹⁰

1. Use 4 µl sample with 300 µl of Direct HDL Cholesterol Reagent 1
2. Equilibrate to 37°C for 5 minutes.
3. Add 100 µl of Direct HDL Cholesterol Reagent 2
4. Equilibrate to 37°C for 5 minutes.
5. Measurement (Absorb. Difference between 700nm & 600 nm)

CALCULATIONS

(A = Absorbance)

$$\frac{A(\text{patient})}{A(\text{standard})} \times \text{Concentration of calibrator} = \text{Concentration of patient (mg/dl)}$$

Example:

$$\frac{A(\text{patient})}{A(\text{calibrator})} = 0.40, A(\text{calibrator}) = 0.32, \text{Concentration of calibrator} = 53 \text{ mg/dl}$$

$$\frac{0.40}{0.32} \times 53 = 66.3 \text{ mg/dl}$$

To convert from conventional units to SI units, multiply the conventional units by 0.02586.

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mg/dl $\times 0.02586 = \text{mmol/L HDL cholesterol}$

LIMITATIONS

1. Anticoagulants containing citrate should not be used.
2. Protect the reagents from direct sunlight.
3. Store the reagents as per instructions.
4. Samples with values greater than 150 mg/dl must be diluted 1:1 with saline and re-assayed. Multiply the result by two.

CALIBRATION

The Direct HDL/DL Cholesterol Calibrator is required for calibration. Calibrate with each bottle change or lot change or if control results are found to be out of range. The value of the Direct HDL/DL calibrator was assigned by procedures traceable to the National Reference System for Cholesterol (NRSCHOL). Refer to Direct HDL/DL Cholesterol calibrator kit package insert for instructions. If control results are found to be out of range, the procedure should be recalibrated.

QUALITY CONTROL

Reliability of test results should be routinely monitored with control materials that reasonably emulate performance of patient specimens. Quality control materials are intended for use only as monitors of accuracy and precision. Any HDL/DL controls would be suitable for use with this assay. The National Cholesterol Education Program (NCEP) Lipid Standardization Panel (LSP) recommends two levels of controls, one in the normal range (35-65 mg/dl) and one near the concentrations for decision making (<35 mg/dl). An acceptable range of HDL cholesterol values should be established for the controls by repeat analysis. The recovery of control values within the appropriate range should be the criteria used in evaluation of future assay performance. Quality control materials are intended for use only as monitors of accuracy and precision. Controls should be run with every working shift in which HDL-C assays are performed. It is recommended that each laboratory establish their frequency of control determination. Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

EXPECTED VALUES

The expected values for serum HDL cholesterol are as follows:¹⁴

- < 40 mg/dl Low
- ≥ 60 mg/dl High

It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

According to the NCEP, HDL values greater than or equal to 35 mg/dl are considered desirable, and values greater than or equal to 60 mg/dl are considered to offer some protection against coronary heart disease. Values below 35 mg/dl are considered to be a significant independent risk factor for coronary heart disease.

PERFORMANCE CHARACTERISTICS

Assay Range: 2 - 150 mg/dl

Accuracy: Accuracy of the Direct HDL Cholesterol Reagent method was verified by comparison to the Designated Comparison Method (ultracentrifugation, chemical precipitation and enzymatic cholesterol analysis which has been

standardized to the Abell-Kendall method)¹⁰ and another manufacturer's lyophilized Direct HDL Cholesterol method to the Designated Comparison Method produced the following results:

Studies comparing the Direct HDL Cholesterol method to Pointe Scientific AutoHDL Cholesterol Reagent Set on Hitachi 717 Analyzer produced the following results:

Method	Direct HDL Cholesterol	AutoHDL Cholesterol
N	52	52
Mean HDL Cholesterol	46	48.3
Range (mg/dl)	22-96.8	21.3-92.1
Regression Analysis	$y=0.93x+0.73$ mg/dl	
Correlation Coefficient	$r=0.916$	

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Precision. Within Day precision for the Direct HDL Cholesterol Reagent was determined on Hitachi 717 following a modification of NCCLS document EP5-T2. Within Day precision studies produced the following results:

	Sample 1	Sample 2
N	25	25
Mean HDL Cholesterol	38.1	83.5
Standard Deviation (mg/dl)	1.54	1.94
Correlation of Variation (%)	4.0	2.3

Day to Day
Day to Day precision for the Direct HDL Cholesterol Reagent was also determined on Hitachi 717 following a modification of NCCLS document EP5-T2. Day to Day precision studies produced the following results:

	Sample 1	Sample 2
N	25	25
Mean HDL Cholesterol	38.2	83.4
Standard Deviation (mg/dl)	1.74	2.2
Correlation of Variation (%)	4.5	2.6

Sensitivity The analytical sensitivity of the Direct HDL Cholesterol Reagent was determined on Hitachi 717 as 2 mg/dl of HDL Cholesterol.

REFERENCES

1. Gotto A.M., Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital practice, 23:Suppl.1,4 (1988).
2. Crouse J.R. et al, Studies of low-density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26:666 (1985).
3. Badimon J.J., Badimon L., Fuster V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 1990; 85:1234-41.
4. Castelli W.P., et al, Cholesterol and other lipids in coronary heart disease, Circulation, 55:767 (1977).
5. Barr, D.P., Russ E.M., Eder H.A., Protein-lipid relationships in human plasma, Am. J. Med. 11:480 (1951).
6. Gordon, T., et al, High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med., 62:707 (1977).
7. Williams, P., Robinson D., Baily A., High-density lipoprotein and coronary risk factor, Lancet, 1:72 (1979).
8. Kannel W.B., Castelli W.P., Gordon, T., Cholesterol in the prediction of atherosclerotic disease, New perspectives based on the Framingham study, Am. Intern. Med., 90:85 (1979).
9. National Institute on Health publication No. 93-3095, September 1993.
10. Warnick G. Russell, Wood Peter D., National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol, Executive Summary Clinical Chemistry, Vol. 41, No. 10 (1995).
11. Grundy, S.M., et al, Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II), JAMA 1993, 269: 23, 3015-3023.
12. National Committee for Clinical Laboratory Standards, National Evaluation .

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High Technology, Inc

PRODUCT NAME:

Liquid Stable Lithium Enzymatic Assay Kit (200 Tests)
Reagent R1: 2 x 20 mL
Reagent R2: 2 x 10 mL
Calibrator: 3 x 3 mL
Control: 2 x 3 mL

Intended use:

Diazyme Lithium Enzymatic Assay kit is for quantitative *in vitro* determination of lithium in human serum.

Storage and Handling:


The enzymatic lithium reagent should be stored at 2-8°C. Shelflife - 1 year.

Product feature:

For decades, lithium carbonate has remained one of the most effective agents for treatment of patients suffering from bipolar disorder (manic depressive psychosis). Recently, studies have also shown that lithium holds promise against Alzheimer's disease. However, lithium has many side effects. Overdosage of lithium can cause acute intoxication, which occurs quite often due to its narrow therapeutic index. For example, serum lithium levels over 1.5 mmol/L (12 hours after a dose) usually indicate a significant risk of intoxication. Therefore, the timely and accurate monitoring of serum levels of lithium after a therapeutic dosage is critical.

Assay Principle

Lithium is determined spectrophotometrically through a kinetic coupled enzyme assay system involving Diazyme's proprietary phosphatase whose activity is sensitive to lithium. Through enzymatic coupling, the phosphatase substrate, adenosine biphosphate (PAP) is converted to hypoxanthine by a series of enzymatic reactions to generate uric acid and hydrogen peroxide. The hydrogen peroxide generated reacts with N-Ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylolindole (EHSPT) and 4-aminopyrene (4-AA) in the presence of peroxidase (POD) to form a quinone dye which has maximal absorbance at 556nm. The rate of the quinone dye formation is inversely proportional to the concentration of lithium in serum samples.

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
Performance

Accuracy:

To demonstrate accuracy, the Diazyme Lithium Enzymatic Assay was tested with individual serum samples compared to a currently marketed Lithium ISE method. The Diazyme Enzymatic Lithium assay method (y) results generated on a Roche Hitachi 917 analyzer were compared with those obtained with a commercially available ISE method (x). A total of sixty two (62) patient samples and 4 controls were tested with values ranging from 0.1 to 3.1 mmol/L. These studies yielded a correlation coefficient of 0.9916 with a linear regression of $y = 1.0344x - 0.0445$. Values above 1.5 mmol/L were supplemented with a stock solution of lithium acetate to demonstrate performance across the entire reportable range.

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
High Technology, Inc


Lithium Reagent Set Литий

Diazyme Lithium Assay - ISE Correlation
(47 unaltered samples, 15 altered serum samples)

Linearity

Diazyme Enzymatic Lithium Assay has a linear range from 0.19-3.0 mmol/L of lithium.

Precision

The precision of the Diazyme Lithium Enzymatic Assay was evaluated according to Clinical and Laboratory Standards.


Within-run Precision	Level 1 1.0 mmol/L	Level 2 2.5 mmol/L
Number of Data Points	40	40
Mean (mM)	0.97	2.50
SD (mM)	0.04	0.03
CV%	4.3	1.2
Total Precision	Level 1 1.0 mmol/L	Level 2 2.5 mmol/L
Number of Data Points	40	40
Mean (mM)	0.97	2.50
SD (mM)	0.05	0.03
CV%	4.8	1.3

Interference

To determine the level of interference from the substances normally present in the serum, Diazyme Liquid Stable Enzymatic Assay was tested with 1.0 mmol/L and 2.5 mmol/L lithium serum samples spiked with various concentrations of substances following Clinical and Laboratory Standards Institute (formerly NCCLS) EP7-A 'Interference Testing in Clinical Chemistry' dose-response guidelines. The following substances normally present in the serum produced less than 10% deviation when tested at levels equal to the concentrations listed.


Interference-Concentration
Ascorbic Acid-5 mM
Bilirubin-45 mg/dL
Bilirubin Conjugate-20 mg/dL (low lithium serum) 45 mg/dL (high lithium serum)
Hemoglobin-500 mg/dL
Triglycerides-1000 mg/dL
Ca²⁺-4 mM
Cu²⁺-0.25 mM
Fe³⁺-0.25 mM
K⁺-10 mM
NH₄⁺-0.5 mM

Diazyme Lithium Assay - ISE Correlation
(47 unaltered samples, 15 altered serum samples)

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Magnesium Reagent Set Магний

INTENDED USE

For the qualitative determination of magnesium in serum

INTRODUCTION

Magnesium is one of the most abundant cations in the body and is essential to many physiological processes. Approximately one-half of the body magnesium is present in the bone. The remainder is found in soft tissues and blood cells with a small amount present in blood. Magnesium is an activator of various enzymes and is also essential for the preservation of the macromolecular structure of DNA, RNA, and ribosomes. *

Little is known about the factors regulating magnesium levels in plasma. It is believed that the parathyroid gland may be involved. Decreased levels have been observed in cases of diabetes, alcoholism, diuretics, hyperthyroidism, malabsorption, hyperparathyroidism, myocardial infarction, congestive heart failure, and liver cirrhosis. Increased serum magnesium levels have been found in renal failure, diabetic acidosis, Addison's disease, and vitamin D intoxication. In terms of accuracy, speed, and convenience, the determination of magnesium by atomic absorption spectrophotometry is the method of choice. However, this method requires expensive instrumentation and uses large sample volumes that limit its use for frequent testing. This procedure is a direct method in which magnesium forms a colored complex with calmagite in a basic solution, where calcium and protein interference is eliminated by EGTA and surfactant. *

PRINCIPLE

Magnesium forms a colored complex with calmagite in alkaline medium to produce a red complex that is measured spectrophotometrically at 530 nm. EGTA serves to complex and prevent calcium interference, and a surfactant eliminates the effect of protein. The color produced is proportional to the magnesium concentration.

REAGENT COMPOSITION

1. Magnesium Buffer Reagent:
2-Ethylaminoethanol 6.0 w/v, potassium cyanide 0.10% w/v, EGTA 1.18 mM. **Caution: Contains Cyanide, avoid all contact and ingestion.**
2. Magnesium Color Reagent:
Calmagite 0.006% w/v, stabilizer 2.0% w/v, surfactant 0.03% w/v. **Caution: DO NOT PIPETTE by mouth.**
3. Magnesium Standard:
2 mEq/L magnesium iodate, tetrahydrate.

WARNINGS AND PRECAUTIONS:

1. For *in vitro* diagnostic use only.
2. Exercise the normal precautions required for the handling of all laboratory reagents. Pipetting by mouth is not recommended for any laboratory reagent.

REAGENT PREPARATION

The working reagent is prepared by mixing ten (10) volumes of color reagent with one (1) volume of buffer reagent in a disposable plastic container. Combine only the volume of reagent necessary to perform the specific number of tests for that day. Working reagent is stable for twenty-four (24) hours at 18-25°C. Disposable plastic containers or acid-washed glass containers are recommended to avoid contamination.

STORAGE AND STABILITY

The reagent kit is stable until the expiration date stated on the label if stored refrigerated between 2 - 8°C. The working reagent is stable for twenty-four (24) hours at 18 - 25°C.

REAGENT DEGRADATION

1. Failure to achieve assayed values on freshly prepared control sera would indicate deterioration.
2. Working reagent becomes visibly turbid.

SPECIMEN COLLECTION

1. Unhemolyzed sera are the recommended sample. Separate from clot as soon as possible.
2. Avoid taking blood from a limb that is simultaneously receiving an infusion.
3. Avoid contamination of blood with tissue fluid.

INTERFERING SUBSTANCES

1. Plasma collected with anticoagulants such as EDTA, citrate and oxalate must not be used.
2. A number of drugs and substances affect the concentration of magnesium. See Young, et al.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Pipetting devices.
2. Test tubes/rack.
3. Timing device.
4. Spectrophotometer capable of reading at 530 nm.

PROCEDURE (AUTOMATED)

Consult the appropriate instrument application guide available from Tecco.

PROCEDURE (MANUAL)

1. Prepare working reagent according to preparation instruction.
2. Label test tubes: "Blank", "Standard", "Control", "Patient", etc. For each sample, dispense 1.0 ml of working reagent to each tube. Add 0.01 ml (10 µl) sample to its respective tube. Mix gently.

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High Technology, Inc

Incubate for five (5) minutes at room temperature.
After incubation, zero spectrophotometer with the reagent blank at 530 nm. (Wavelength range: 500-550 nm).
Read and record absorbances of samples.

TO

* TC - MULTI PURPOSE CALIBRATOR MAY BE USED REPLACE STANDARD.

NOTE: Final color is stable for at least thirty-minutes at room temperature.

CALCULATIONS

Abs = Absorbance

$\frac{\text{Abs. of Unknown}}{\text{Abs. of Standard}} \times \text{concentration of standard} = \text{concentration of unknown in mEq/L}$

Example: Abs. of unknown = 0.098
Abs. of standard = 0.113
Concentration of standard = 2.0 mEq/L
Then $0.098 \times 2 \text{ mEq/L} = 1.73 \text{ mEq/L}$
0.113

Quality Control

It is recommended that controls be included in each set of assays. Commercially available control material with established magnesium values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction, or procedural errors.

Expected Values

Adults 1.3-2.5 mEq/L

Linearity: 4.0 mEq/L (4.86 mg/dl)

Comparison: Studies between the present method and a similar calmagite method yielded a coefficient correlation of 0.96 with a regression equation of $y = 0.98x + 0.09$ sample values ranged from 1.5 mEq/L to 2.8 mEq/L.

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Copper Reagent Set Медь

Product Summary

Stability :	Until Expiry at 18-25°C
Linear Range :	3-85 $\mu\text{mol/L}$ (19-540 $\mu\text{g/dL}$)
Specimen Type :	Serum
Method :	Endpoint
Reagent Preparation :	Mix equal quantities of Reagent 1 and Reagent 2.

INTENDED USE

This reagent is intended for the *in vitro* quantitative determination of Copper in serum on automated clinical chemistry analyzers.

CLINICAL SIGNIFICANCE

Copper is an essential trace element found predominantly bound to the copper transport protein ceruloplasmin, while a very small proportion is complexed with albumin and other metalloproteins.

The most significant clinical application of copper determination is in the diagnosis of Wilson's disease. This is associated with a decrease in the synthesis of ceruloplasmin, which results in low serum copper levels. A second disorder of copper metabolism is Menkes' syndrome or kinky hair syndrome, an X-linked genetic defect in copper absorption. Low serum copper levels have also been observed in a number of hypoproteinemias, while increased levels are found in a number of acute and chronic diseases such as leukemia, hemochromatosis and biliary cirrhosis.

METHODOLOGY

Copper bound to ceruloplasmin is released by the reducing agent Guanidine hydrochloride in an acidic medium. 2-(5-bromo-2-pyridylazo)-5-(N-propyl-N-sulfo)pyrimidine (5-Br-PSAA) reacts with the free copper to form a stable colored complex. The intensity of this color is proportional to the copper concentration in the sample and is measured photometrically at 580 nm.

REAGENT COMPOSITION

Active Ingredients Concentration

Reagent 1:

Acetate buffer, pH 4.2 0.4 mol/L
Guanidine HCl 5 mol/L
Catalyst
Oxidisable agent

Reagent 2:

5-Br-PSAA, sodium salt 0.1 mmol/L

WARNING: Avoid ingestion and contact with skin, mouth and eyes. The toxicity of this reagent has not been determined. Flush with plenty of water when disposing. For further information, consult the Copper Reagent Material Safety Data Sheet.

R22 Harmful if swallowed.

R36/38 Irritating to eyes and skin.

S23 Do not breathe vapour.

REAGENT PREPARATION

Prepare the working reagent by mixing equal quantities of Reagent 1 and Reagent 2.

2. Refer to specific instrument applications, available upon request from the Technical support group.

STABILITY AND STORAGE

Prior to Use:

When stored between 18-25°C, the reagents are stable until the expiration date stated on the bottle and kit box labels.

Working Reagent:

Working reagent is stable for at least 5 days at 2-8°C.

Indications of Reagent Deterioration:

- Turbidity;
- Reagent absorbance > 0.2 AU at 580 nm; and/or
- Failure to recover control values within the assigned ranges.

SPECIMEN COLLECTION AND HANDLING

Serum: Use non-hemolyzed serum.

Storage: Serum copper samples are stable for at least 2 days at room temperature (18-25°C) or for 1 week at 2-8°C.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyser capable of maintaining constant temperature

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Copper Reagent Set Медь

(37°C)

and measuring absorbance between 550 and 590nm.

- It required, pipettes for accurately dispensing measured volumes.
- Analyser specific consumables, e.g. sample cups.
- Assayed Normal and Abnormal control material.
- A suitable aqueous Copper standard.

ASSAY PROCEDURE

The following system parameters are recommended. Individual instrument applications are available upon request from the Technical Support Group.

SYSTEM PARAMETERS

Temperature 37°C
Primary Wavelength 580 nm (550 - 590nm)
Assay Type End Point
Direction Increase
Sample:Reagent ratio 1:10
e.g. Sample vol 0.1 mL
Combined Reagent vol 1.0 mL
Incubation Time 300 seconds
Reagent Blank Limits Low 0.10 AU
(580nm, 1cm lightpath) High 0.20 AU
Linearity 3 - 85 $\mu\text{mol/L}$ (19 - 540 $\mu\text{g/dL}$)
Sensitivity 4.1 ΔmAbs per $\mu\text{mol/L}$
(550 nm, 1 cm lightpath) 0.7 ΔmAbs per $\mu\text{g/dL}$

CALCULATIONS

Results are calculated, usually automatically by the instrument, as follows:-
Copper = Abs of Unknown/ Abs of Calibrator x Calibrator Value

Example:

Absorbance of calibrator = 0.081
Absorbance of unknown = 0.136
Value of calibrator = 15.73 $\mu\text{mol/L}$ (100 $\mu\text{g/dL}$)
Copper = $0.136 / 0.081 \times 15.73 = 26.4 \mu\text{mol/L}$
Copper = $0.136 / 0.081 \times 100 = 168 \mu\text{g/dL}$

NOTES

1. The reagent and sample volumes may be altered proportionally to accommodate different spectrophotometer requirements.
2. Specimens with copper values above 85 $\mu\text{mol/L}$ (540 $\mu\text{g/dL}$) should be diluted with isotonic saline and reassayed. Multiply results by the dilution factor.
3. The color reaction is stable for at least 10 minutes at 37°C.
4. Unit Conversion: $\mu\text{mol/L} \times 6.355 = \mu\text{g/dL}$

CALIBRATION

Calibration is required. A suitable aqueous copper standard is recommended. For calibration frequency on automated instruments, refer to the instrument manufacturer's specifications. However, calibration stability is contingent upon optimum instrument performance and the use of reagents which have been stored as recommended in the stability and storage section of this package insert. Recalibration is recommended at anytime if one of the following events occurs:-

- The Lot number of reagent changes
- Preventative maintenance is performed or a critical component is replaced
- Control values have shifted or are out of range and a new vial of control does not rectify the problem.

QUALITY CONTROL

To ensure adequate quality control, normal and abnormal control with assayed values should be run as unknown samples:-

- At least every eight hours
 - When new bottles of reagent are used
 - After preventative maintenance is performed or a critical component is replaced
 - With every calibration
- Control results falling outside the upper or lower limits of the established ranges indicate the assay may be out of control. The following corrective actions are recommended in such situations:-
- Repeat the same controls.
 - If repeated control results are outside the limits, prepare fresh control serum and repeat the test.
 - If results are still out of control, recalibrate with fresh standard, then repeat the test.
 - If results are still out of control, perform a calibration with fresh reagent, then

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High Technology, Inc

repeat the test.

• If results remain out of control contact Technical Services or your local distributor.

LIMITATIONS

1. Studies to determine the level of interference from iron and zinc were carried out. The following results were obtained:
Iron: No interference from iron up to 100 $\mu\text{mol/L}$.
Zinc: No interference from zinc up to 100 $\mu\text{mol/L}$.
2. Avoid lipaemic, haemolysed and icteric samples.
3. Young DS has published a comprehensive list of drugs and substances which may interfere with this assay.

EXPECTED VALUES

Males 11.0 - 22.0 $\mu\text{mol/L}$ (70 - 140 $\mu\text{g/dL}$)

Females 12.6 - 24.4 $\mu\text{mol/L}$ (80 - 155 $\mu\text{g/dL}$)

The quoted values are representative of the expected range for this method and should serve as a guide only. It is recommended that each laboratory verify this range or derive a reference interval for the population that it serves.

PERFORMANCE DATA

The following data was obtained using the Copper reagent on a well maintained automated clinical chemistry analyser. Users should establish product performance on their specific analyser used.

IMPRECISION

Imprecision was evaluated using two levels of commercial control and following the NCCLS EP5-T procedure.

Within run:	LEVEL I	LEVEL II
Number of data points	80	80
Mean ($\mu\text{mol/L}$ / $\mu\text{g/dL}$)	16.25 / 103.27	31.50 / 200.00
SD ($\mu\text{mol/L}$ / $\mu\text{g/dL}$)	0.46 / 2.92	0.56 / 3.56
CV (%)	2.83	1.78

Total:	LEVEL I	LEVEL II
Number of data points	80	80
Mean ($\mu\text{mol/L}$ / $\mu\text{g/dL}$)	16.25 / 103.27	31.50 / 200.00
SD ($\mu\text{mol/L}$ / $\mu\text{g/dL}$)	0.80 / 5.08	1.03 / 6.55
CV (%)	4.92	3.27

ACCURACY

Comparison studies were also carried out using a similar commercially available Copper reagent as a reference. Serum samples were assayed in parallel and the results compared by least squares regression. The following statistics were obtained:

Number of sample pairs 50
Range of sample results 5 - 35 $\mu\text{mol/L}$ (32 - 222 $\mu\text{g/dL}$)
Slope 0.95
Intercept -0.43 $\mu\text{mol/L}$ (2.7 $\mu\text{g/dL}$)
Correlation coefficient 0.997
Comparison studies were also carried out using AAS as a reference. The following statistics were obtained:
Number of sample pairs 60
Range of sample results 3.5 - 35.0 $\mu\text{mol/L}$ (22 - 222 $\mu\text{g/dL}$)
Slope 0.96
Intercept 1.15 $\mu\text{mol/L}$ (7.3 $\mu\text{g/dL}$)
Correlation coefficient 0.97

LINEARITY

When run as recommended, the assay is linear between 3 and 85 $\mu\text{mol/L}$ (19 and 540 $\mu\text{g/dL}$).

Linearity on various automated instruments may vary from this value. The user should consult the specific instrument application for the instrument specific linearity claim.

SENSITIVITY

When run as recommended, the sensitivity of this assay is 4.1 ΔmAbs per $\mu\text{mol/L}$ or approximately 0.7 ΔmAbs per $\mu\text{g/dL}$ (1cm light path, 550 nm).

REFERENCES

1. Tietz NW (Ed). "Textbook of Clinical Chemistry" WB Saunders, 1986; 929-933.
2. Makino T., Clin Chem Acta, 1989; 185:7-16.
3. Horiguchi D et al. Anal. Chim. Acta, 1983; 151:457-463.
4. Young DS, Effects of Drugs on Clinical Laboratory Tests, Third Edition, 1990.
5. Wiedel M et al. Creation and Verification of Reference Intervals. Laboratory Medicine 1995; 26:593-7.
6. National Committee for Clinical Chemistry Standards. User evaluation of Precision Performance of Clinical Laboratory Devices NCCLS, 1984, NCCLS Publication

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Myoglobin (Mb) Reagent Set Миоглобин

Intended use:

To measure the level of myoglobin in human serum or plasma by latex immunoturbidimetric method.

The principle of the method:

Latex reagent is a suspension of polystyrene latex particles coated by anti-human myoglobin antibodies. Mixing human serum or plasma containing myoglobin with a latex reagent results in agglutination of the particles accompanied by increase in optical density that can be quantified by absorbance at 570 nm. Concentration of myoglobin in the sample is determined using a standard curve constructed by measuring the optical density of standards.

Reagents:

Reaction Buffer: glycine buffer

Latex reagent: 0.12% suspension of latex particles coated with antibodies to human myoglobin.

Reagent preparation:

All reagents are ready-to-use

Storage:

+2 +4 °C

Cautions

1. Only for in vitro diagnostics.
2. Not for ingestion.
3. Avoid contact with reagent. In the case of contact, wash with water.
4. Reagent contains sodium azide as preservative. Sodium azide may form explosive compounds when contact with the metals.

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Intended Use

For the quantitative determination of low levels of albumin in urine by immunoturbidimetric assay. For *in vitro* diagnostic use.

Introduction

A small amount of protein is excreted daily into the urine of healthy individuals. The excreted proteins are mucoproteins, most of which are filtered out of the uriniferous tubules and the glomeruli. Albumin, a protein of molecular weight of 50,000, is not easily filtered out and is excreted into the urine (microalbuminuria) 1,2. This makes albumin excretion into the urine a useful indicator of early glomerular disease. Microalbuminuria is a condition characterized by increased urinary excretion of albumin in the absence of overt nephropathy 3,4. Microalbumin has been reported in several studies to predict development of diabetic nephropathy and its mortality risk in diabetic patients. 2,5,6,7,8,9. Because microalbuminuria may be reversible if diabetes is well controlled, the early detection of microalbumin may be very beneficial in treatment programs for diabetes. Albumin in urine has been measured by a variety of methods. Many of these lack the sensitivity to be used for microalbumin measurements. Radioimmunoassay (RIA) and immunoturbidimetric assay provide the sensitivity required 10. This Microalbumin uses an immunoturbidimetric format.

Principle

When a sample is mixed with anti-human albumin goat antiserum, agglutination is caused by the antigen-antibody reaction. The turbidity is measured at 340nm and 700nm and albumin in the sample is quantitatively determined.

Characteristics

1. No reagent preparation such as dissolving, mixing, or dilution.
2. No sample dilution up to 30 mg/dl (300 µg/ml).
3. Long reagent stability after opening (1 month).
4. Good correlation with RIA.
5. No interference from drugs.
6. No prozone effect in ordinary measuring range.
7. Compatible with most clinical chemistry auto-analyzers.

Reagents (Liquid stable)

- R1: Buffer Reagent, pH 7.6
100mM Tris (hydroxymethyl) aminomethane
R2: Antiserum Reagent, pH 7.6
20% Anti-human albumin, goat antiserum
100mM Tris (hydroxymethyl) aminomethane

Warnings and Precautions

1. FOR IN VITRO DIAGNOSTIC USE.
2. Not to be used internally in humans or animals.
3. Normal precautions exercised in handling laboratory reagents should be followed.
4. Do not mix or use reagents from one test kit with those from a different lot number.
5. Do not use reagents past their expiration date stated on each reagent container label.
6. Do not pipette by mouth. Avoid ingestion and contact with skin.

Reagents in this kit contain sodium azide as a preservative. Sodium azide may form explosive compounds in metal drain lines. When disposing of reagents through plumbing fixtures, flush with copious amounts of water.

Reagent Preparation

Reagents are ready to use and do not require reconstitution.

Storage and Handling

All reagents should be stored refrigerated (2-10°C). Return all reagents to 2-10°C promptly after use. Unopened reagents can be used up to the expiration date on the package and bottle labels.

Reagent Stability

Discard reagents if they become contaminated. Evidence of cloudiness or particulate material in solution is cause to discard. If the absorbance of isotonic saline is greater than 0.1 or if the absorbance of the 5.0 mg/dl calibrator is smaller than 0.15, the reagents should not be used. Opened reagents can be used for 1 month if stored at 2-10°C.

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The specimen should be a fresh or a 24-hour urine. Urine specimens should

Be stored refrigerated (2-10°C). The specimens may be stored refrigerated up to two weeks or frozen at -70°C for at least 5 months. 11

Automated Analyzer Application

Suitable for two-reagent automated analyzers that use a two-point calibration method. Measurements of absorbance are to be made with a spectrophotometer able to accurately read absorbance at 340 and 700nm. Refer to the instrument manual from the manufacturer regarding the following:

- Use or function
- Installation procedures and requirements
- Principles of operation
- Performance characteristics and specifications
- Operation instructions
- Calibration procedures including materials and/or equipment to be used
- Operational precautions and limitations
- Hazards
- Service and maintenance

Materials Provided

- Reagent 1 (R1) Buffer Reagent 4x20ml
- Reagent 2 (R2) Antiserum Reagent 2x10ml

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- Calibrators: Microalbumin Multi-Calibrator Set, 6 Calibrators, Approx. values: 0, 0.5, 1.0, 5.0, 10.0, 30.0 mg/dl (For actual values see vial labels)
- Spectrophotometer: capable of accurate absorbance reading at 340 and 700nm with appropriate cuvettes.
- Pipettes: capable of accurately dispensing the required volumes
- Test Tubes: glass or plastic
- Water Bath: capable of maintaining 37°C

Procedure (Automated)

Refer to specific instrument application instructions.

Calibration Curve

For albumin concentrations less than 10 mg/dl, a two-point calibration curve can be made using a saline blank (0mg/dl) and an albumin standard. When a sample with a higher albumin concentration (>10mg/dl) is assayed, it is recommended that a multi-point calibration curve using the microalbumin multi-calibrator set be made. It is recommended that a calibration curve be made each day.

Quality Control

A quality control program is recommended for all clinical testing laboratories. It is recommended that control urines, both normal and abnormal, be run with each batch of samples to monitor the procedure.

The values obtained for controls should fall within the manufacturers specified range. A laboratory may establish its own control urine by assaying the urine a sufficient number of times to generate a valid mean and acceptable range.

Calculations

Albumin levels are determined using the prepared calibration curve.

Limitations

The measurable range for this albumin test kit is between 0.5 mg/dl and 30 mg/dl. If albumin concentrations are greater than 30 mg/dl, dilute 1 part sample with 4 parts isotonic saline including 0.5% Tween 20 and re-assay. Multiply the result by 5 to compensate for the dilution.

Performance

Sensitivity

When saline blank is used as a sample, the absorbance is below 0.05. When a calibrator containing 5 mg/dl of human albumin is assayed, the absorbance (after subtracting the absorbance value for the saline blank) is within range of 0.131-0.525.

Specificity

When a urine sample with a known value is assayed, it is within $\pm 10\%$.

Precision:

When a sample containing 5 mg/dl of human albumin is

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repeatedly assayed 20 times, the absorbance C.V. is less than 5%.

Precision Assay:

Within Run (N=20)			Run to Run (N=20)		
Mean	S.D.	C.V.%	Mean	S.D.	C.V.%
0.96	0.05	5.2	0.97	0.11	11.3
5.47	0.12	2.2	5.26	0.22	4.2

Assay Range

0.5-30 mg/dl or 5-300 ug/ml (multi-point calibration)

0.5-10 mg/dl or 5-100 ug/ml (two-point calibration)

Correlation:

A comparison of this Microalbumin and a Kamiya Microalbumin Test Kit was performed on a Hitachi 717 automated analyzer. The test results provided the following data. All values are expressed in mg/dl of albumin.

$$y = 1.0738x + 0.01555$$

$$r = 0.995, (n=74, \text{ range} = 0.9-20.7)$$

x = Kamiya Test Kit y = this Microalbumin test kit

x = min = 0.900 y min = 1.000

x max = 20.7 y max = 22.7

x mean = 5.86 y mean = 6.33

Interferences

Ascorbic acid: No interference up to 200 mg/dl (Less than 5%)

Glucose: No interference up to 3.0 g/dl (Less than 15%)

Uric Acid: No interference up to 100 mg/dl (Less Than 5%)

Creatinine: No interference up to 300 mg/dl (Less than 8%)

Creatine: No interference up to 100 mg/dl (Less than 5%)

Calcium: No interference up to 30.0 mM (Less than 8%)

NaCl: No interference up to 900 mg/dl (Less than 15%)

Mg: No interference up to 30 mM (Less than 5%)

KCl: No interference up to 300 mg/dl (Less than 5%)

Urea: No interference up to 3.0 g/dl (Less than 8%)

Expected Values

The expected value for Microalbumin is 30-300 mg/24 hours. Each laboratory should establish its own expected values using this kit.

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Intended Use

For the quantitative determination of total protein in urine.

Clinical Significance

The presence of protein in urine is a very sensitive indicator of renal disorders. There are four ways by which increased amounts of protein can occur: increased glomerular permeability, defective tubular re-absorption, increased plasma concentration of an abnormal, low molecular weight protein, and abnormal secretion of protein into the urinary tract. Albuminuria, increased amounts of albumin in urine, has been recognized as an early indicator of renal damage in diabetes that can be reversed if detected and treated early. 2

Method History

Various methods have been described for the determination of protein concentrations in biological fluids. These methods are based on colorimetric, turbidimetric, electrophoretic, or immunologic procedures. 3,4 The dye binding methods are characterized as having good precision and sensitivity. The Coomassie Blue method is very sensitive, but the reagents stain glassware and plastic.

This method is based on the procedure developed by Fujita and Watanabe. 7 It is a sensitive dye binding, colorimetric method employing Pyrogallol Red. The method seldom stains cuvettes or plastic tubing, and can be automated.

Principle

Pyrogallol Red is combined with molybdenum acid at a low pH. When the complex is combined with protein, a blue-purple color is formed. The increase in absorbance at 600 nm is directly proportional to the protein concentration in the sample.

Reagents

1. MICROPROTEIN REAGENT: Contains Buffer, Pyrogallol Red 0.067 mmol/L, sodium molybdate stabilizer 0.153 mmol/L, surfactants, and preservative.

2. PROTEIN STANDARD SOLUTION: Contains Albumin 50 mg/dl in saline with sodium azide 0.05% as a preservative.

Reagent Preparation

Microprotein Reagent and Protein Standard Solution are provided ready to use.

Reagent Storage

Store Microprotein Reagent and Protein Standard refrigerated (2-8°C). Reagent and standard are stable until the expiration date shown on the label.

Precautions

1. Microprotein Reagent is for *in vitro* diagnostic use only.
2. Normal precautions exercised in handling laboratory reagents should be followed.
3. Protein Standard contains sodium azide. Do not ingest. May react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Reagent / Standard Deterioration

The reagent and standard should be clear. Do not use if turbid. Failure to achieve assayed control values may indicate deterioration of reagent and/or standard. Do not use if the reagent absorbance at 600nm is less than 0.100.

Specimen Collection and Storage

It is recommended that specimen collection be carried out in accordance with NCCLS Document M29-T2. Specimens containing visible particulate matter should be clarified by centrifugation prior to testing.

URINE: Tests are performed on 24-hour samples. The urine should not be collected during periods of exercise because of its effect on albumin concentration.

Protein determinations should be performed with fresh specimens. If test cannot be performed with fresh urine, specimens may be stored at -20°C for up to one year. 8 NOTE: Hemoglobin can increase recovered total protein values. DO NOT USE samples containing blood.

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It is recommended not to use urine specimens with added preservatives since

some added preservatives such as HCL and benzoic acid have been shown to interfere in the protein assay, giving false low results. 7 Bilirubin to a level of 20 mg/dl and Ascorbic acid to a level of 3.0 mg/dl have been found not to interfere with the assay. (Less than 3.0% deviation for samples in the range of 130.0-132.0 mg/dl and less than 17% for samples measured in the range of 9.2-12.5 mg/dl) Hemoglobin being a protein, will increase recovered total protein values. Ph variation was found to have no effect on the total protein determination. The effects of specific gravity variation were not evaluated. Some drugs and medications may interfere, see Fujita

Materials Provided

1. Microprotein Reagent
2. Protein Standard Solution

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1. Accurate Pipetting devices
2. Test Tubes/Rack
3. Water Bath or Heat Block (37°C)
4. Timer
5. Spectrophotometer able to read at 600nm

Procedure (Automated)

Application procedures are available for various discrete automated instruments. Please contact Technical Service.

Procedure (Manual)

1. Label test tubes "Blank", "Standard", "Control", "Samples", etc.
2. Pipette 1.0ml of Microprotein reagent to each tube
3. Allow the tubes to warm to 37°C
4. Pipette 0.02ml (20µl) deionized water, standard, controls, and samples to the appropriately labeled tubes
5. Allow tubes to incubate at 37°C for 5 minutes
6. After 5 minutes, set the spectrophotometer to 600nm and zero the instrument with the BLANK tube
7. Read and record the absorbance (Abs) of the standard, controls, and samples
8. To determine Microprotein concentration of the samples, refer to the "Calculations" section.

Calibration

Use an aqueous protein standard (traceable to NIST SRM 927c).

Quality Control

Standard practice for Quality Control should be applied to this procedure. Commercially available Urine and CSF controls (2 levels) should be used to monitor the daily acceptable variations. A satisfactory level of performance is achieved when the analyte values obtained are within the acceptable range established by the laboratory.

Calculations

Protein values are expressed as mg/dl.

Protein (mg/dl) = Abs Unk x Conc. of Std.

Abs Std.

Where:

Abs Unk = The absorbance of the unknown sample

Abs Std. = The absorbance of the standard

Conc. of Std. = Concentration of standard (50 mg/dl)

Example:

Abs Unk = 0.085

Abs Std = 0.195

Conc. of Std. = 50 mg/dl

Protein (mg/dl) = $0.085 \times 50 = 21.8$ mg/dl

0.195

To determine the 24-Hour Urinary Protein, measure the 24-hour urine total volume in ml (TV) and assay the urine protein content (mg/dl). Calculate the 24-Hour Urinary Protein using the following formula:

Protein (mg/day) = Protein (mg/dl) x TV

100

Where: TV = 24-hr. urine total volume in ml

100 = converts ml/day to dl/day

S.I. Units

To convert the results into S.I. units, multiply the microprotein concentration (mg/dl) by 0.0100. For example, microprotein concentration = 21.8 mg/dl x 0.0100 = 0.218 g/L

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NOTE: Each laboratory should confirm the validity of the interval ranges listed for the population it serves.

Performance

Assay Range: The Microprotein procedure has an assay range of 2.0 - 250 mg/dl. Samples that exceed 250 mg/dl should be diluted with an equal volume of isotonic saline and re-assayed. Multiply the result by 2 to compensate for the dilution.

Linearity was performed using eleven known concentrations ranging from 1.23 to 259.43 mg/dl. Study was performed on a Roche Hitachi 917 analyzer.

Sensitivity

Based on an instrument absorbance resolution of 0.001 this procedure has a sensitivity of 0.250 mg/dl. The Limit of Detection was found to be 2.0 mg/dl when the absorbance was read bichromatically at 600 / 700 nm (Hitachi 917 chemistry analyzer).

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INTENDED USE

For the quantitative determination of uric acid in serum.


INTRODUCTION

Uric acid is the end product of purine metabolism. Nearly half of the total uric acid is eliminated and replaced each day by way of urinary excretion and through microbial degradation in the intestinal tract. Increased uric acid levels are commonly associated with both nitrogen retention and urea, creatine, and other non-protein constituents. The quantitation of uric acid is an aid in the diagnosis of gout, decreased renal function, myeloproliferative disorders, and other conditions in which the cause for the hyper-uricemia is not well known.¹

Uric acid is most commonly determined by a phosphotungstate method² and iron reduction method.³ Due to serum interferences, the enzyme uricase has been widely used instead. Uricase is more specific for uric acid since uricase acts only on uric acid.

PRINCIPLE

The enzymatic reaction sequence employed in the assay of uric acid is as follows:

Uric acid is converted by uricase into allantoin and hydrogen peroxide. The hydrogen peroxide initiates the coupling of 4-aminoantipyrine to 3,5-dichloro-2-hydroxybenzene sulfonic acid (DHBS) to form the chromogen which is measured at 520nm and which is proportional to the amount of hydrogen peroxide generated from uric acid.

REAGENTS

1. Uric acid reagent: 4-Aminoantipyrine 4mM, 3,5 Dichloro-2- hydroxybenzenesulfonate 2mM, Stabilizer and Surfactant, buffer pH 7.5.
2. Uric acid standard (5 mg/dl).

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use.
- CAUTION:** The reagents may be hazardous. Handle in accordance with good laboratory procedures, which dictate avoiding ingestion, and eye or skin contact.
2. Serum specimens should be considered infectious and handled appropriately.

STORAGE AND STABILITY

The reagent set is stored refrigerated (2 - 8°C). DO NOT FREEZE. Bring reagent to room temperature before use.

REAGENT DETERIORATION

- The reagent should be discarded if:
1. Turbidity has occurred, turbidity may be a sign of contamination.
 2. There is evidence of discoloration. A slight pink color is normal.

SPECIMEN COLLECTION

1. Test specimen should be serum and free from hemolysis.
2. Bacterial contamination should be avoided to preserve the loss of uric acid.
3. Uric acid in serum is stable for three (3) days at 2 - 8°C and up to six (6) months when frozen.

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Каталожный №
HTI-U582-125
HTI-U582-500
HTI-U582-1000
HTI-U7582-125
HTI-U7582-500
HTI-U7582-1000

INTERFERING SUBSTANCES

1. Bilirubin and ascorbic acid can result in falsely depressed uric acid levels.
2. Lipemic samples may cause falsely elevated uric acid levels.
3. Collection tubes containing formaldehyde as a preservative must be avoided.
4. For a comprehensive review of drug interferences refer to Young et al.7

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1. Pipette devices
2. Test tubes/rack
3. Timing device
4. Heating block (37°C)
5. Spectrophotometer capable of reading at 520 nm

GENERAL INSTRUCTIONS

The reagent for uric acid is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

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MANUAL PROCEDURE

1. Label test tubes, "reagent blank", "standard", "control", "unknown", etc.
2. Pipette 1.0 ml of working reagent into all tubes.
3. Pre-warm all tubes at 37°C for three (3) minutes.
4. Add 0.025 ml (25 µl) of sample to respective tubes and mix.
5. Incubate all tubes at 37°C for ten (10) minutes.
6. After incubation, zero the spectrophotometer with the reagent blank at 520 nm and read/record the absorbance of all tubes. (Wavelength range: 500 – 550 nm).
7. Repeat procedure for each sample.

* TC MULTI-PURPOSE CALIBRATOR MAY BE USED TO

REPLACE STANDARD

ALTERNATE VOLUMES

If the spectrophotometer being used requires a final volume greater than 1.0 ml for accurate reading, use 0.05 ml (50 µl) of sample to 3.0 ml of Reagent. Perform the test as described above.

PROCEDURAL LIMITATIONS

The reagent is linear to 25 mg/dl uric acid. Samples with values exceeding 25 mg/dl should be diluted 1:1 with saline, re-assayed and the results multiplied by two (2). Lipemic samples will give falsely elevated results and a serum blank must be run.
 Serum blank: Add 0.025 ml (25 µl) of sample to 1.0 ml water. Zero the spectrophotometer with water. Read and record absorbance and subtract reading from test absorbance.

CALCULATIONS (RATIOMETRIC)

A = Absorbance

$\frac{A_{\text{unknown}}}{A_{\text{standard}}} \times \text{concentration} = \text{value for unknown (mg/dl)}$

A standard of standard

Example:

If the unknown A = 0.170, standard A = 0.180, concentration standard = 5 mg/dl, then:
 $\frac{0.170}{0.180} \times 5 = 4.7 \text{ mg/dl}$

SI UNITS (mmol/L). Multiply the result (mg/dl) by 10 to convert

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dl to liter and divide by 168 (the molecular weight of uric acid).
mg/dl x 10 = mmol/L, mg/dl x 0.0595 = mmol/L.

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established uric acid values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedural errors.

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1. Linearity: 25 mg/dl

2. Sensitivity

Based on an instrument resolution of 0.001 absorbance, the present procedure has a sensitivity of 0.03mg/dl.

3. Comparison: A comparison with another commercial enzymatic uric acid procedure yielded a correlation coefficient of 1.00 with a regression equation of $y = 1.02x - 0.22$.

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INTENDED USE

For the colorimetric determination of sodium in human serum and plasma on TC Matrix analyzers.

PRINCIPLE(1)

Sodium is determined enzymatically via sodium dependent β -galactosidase activity with ONPG as substrate. The absorbance at 405 nm of the product O-nitrophenyl is proportional to the sodium concentration.

ONPG = o-nitrophenyl- β -D-galactopyranoside

SAMPLE

Serum, plasma treated with lithium heparinate

REAGENT CONTENTS:

R1 Buffer 2 x 40 ml

1a Lyophilized Powder 2 x 40 ml

R2 Diluent 1 x 39 ml

2a Lyophilized Powder 1 x 39 ml

Instruction Insert

REAGENT PREPARATION

No preparation is required

REAGENT COMPOSITION

Contents Initial Concentration of Solutions

R1. Buffer/Enzymes

Tris buffer 450 mmole/L, pH 8.0

Cryptand 5.4 mmole/L

β -galactosidase 20.8 U/ml

R2. Diluent/Substrate

Tris buffer 10.0 mmole/L, pH 6.0

O-nitrophenyl galactoside 5.5 mmole/L

Low Calibrator Lot Specific value stated on vial

High Calibrator Lot Specific value stated on vial

SAFETY PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents. Solution 1 contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes, or if ingested, seek immediate medical attention.

Sodium Azide reacts with lead and copper plumbing to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide. Health and Safety Data Sheets are available on request.

The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.

STABILITY AND PREPARATION OF REAGENTS

R1. Buffer/Enzymes

Transfer the contents of one bottle 1a (Lyophilized Powder) to one bottle R1 (Buffer) and dissolve by swirling gently ensuring that the entire contents are completely dissolved. Avoid the formation of foam. Stable for 2 weeks at 2 to 8°C or 5 days at 15 to 25°C.

R2. Diluent/Substrate

Transfer the contents of one bottle 2a (Lyophilized Powder) to one bottle R2 (Diluent) and dissolve by swirling gently ensuring that the entire contents are completely dissolved. Stable for 4 weeks at 2 to 8°C or 2 weeks at 15 to 25°C.

ASSAY PROCEDURE

1. Reconstitute R1 and R2 reagents as described in Reagent Preparation section and keep the reconstituted R1 and R2 reagents on ice bath.

2. In a cuvette, mix 25 μ L of Reagent R1 and 10 μ L of serum sample and incubate at 37°C for 5 minutes and then add 100 μ L of Reagent 2.

3. Read absorbance (405nm) at 3 minutes after addition of

Reagent 2 as A1. incubate for 2 more minutes and read the absorbance as A2.

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
PROCEDURE NOTES

When Sodium and Potassium are requested together Sodium is assayed immediately before Potassium.

PROCEDURE

TEST NAME: NA R1: 250
TEST NO: R2: 100
FULL NAME: Sodium SAMPLE VOLUME: 10
REFERENCE NO: R1 BLANK /
ANALY TYPE: Fixed-time MIX REAG BLANK
PRI WAVE: 405 nm CONCENTRATION /
SECON WAVE: / LINEARITY LIMIT: 80-180
TREND: Ascending SUBSTRATE LIMIT: /
REACT. TIME: 1 - 15 FACTOR: /
INCUBATE TIME: 10 PROZONE CHECK: /
UNIT: mmol/L, Q1: / Q2: / Q3: / Q4: /
PRECISION: Integer PC: / ABS: /
CALIBRATION

The use of Low and High Calibrators is recommended for calibration. A 2-point calibration is recommended every day, with change of reagent lot/bottle or as indicated by quality control procedures.

Ensure that, on the Calibration Checks screen, the following are selected for this test:

Sampling Method for Calibrators

- Duplicate
- Reagent Blank measurement
- Enable Reagent Blank

• Daily
• Reagent Blank System Water
QUALITY CONTROL

Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guideline concerning the running of external quality control. Normal Level and High Level of serum controls are recommended for daily quality control. The two levels of controls should be assayed at least once a day. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error the following steps should be taken:

1. Check instrument settings and light source.
2. Check cleanliness of all equipment in use.
3. Check water, contaminants in bacterial growth may contribute to inaccurate results.
4. Check reaction temperature.
5. Check expiry date of kit and contents.

NORMAL VALUES(2)

136 - 146 mmole/L (313 - 336 mg/dL) It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

LINEARITY

This method is linear between sodium concentrations of 80 and 180 mmole/L (184 and 414 mg/dL).

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INTENDED USE

For the colorimetric determination of sodium in human serum and plasma.

INTRODUCTION

Sodium is the major cation of extracellular fluid. It plays a central role in the maintenance of the normal distribution of water and the osmotic pressure in the various fluid compartments. The main source of body sodium is sodium chloride contained in ingested foods. Only about one-third of the total body's sodium is contained in the skeleton since most of it is contained in the extracellular body fluids.^{1,2}

Hyponatremia (low serum sodium level) is found in a variety of conditions including the following: severe polyuria, metabolic acidosis, Addison's disease, diarrhea, and renal tubular disease. Hypernatremia (increased serum sodium level) is found in the following conditions: hyperadrenalemia, severe dehydration, diabetic coma after therapy with insulin, excess treatment with sodium salts.^{1,2}

PRINCIPLE

The present method is based on modifications of those first described by Maruna and Trinder in which sodium is precipitated as the triple salt, sodium magnesium uranyl acetate, with the excess uranium then being reacted with ferrocyanide, producing a chromophore whose absorbance varies inversely as the concentration of sodium in the test specimen.

REAGENT COMPOSITION

1. **Filtrate Reagent:** Uranyl Acetate 2.1 mM and Magnesium Acetate 20 mM in ethyl alcohol.
2. **Acid Reagent:** A diluted acetic acid.
3. **Sodium Color Reagent:** Potassium Ferrocyanide, non-reactive stabilizers, and fillers.
4. **Sodium Standard:** Sodium Chloride solution: 150 mEq/L of sodium

REAGENT DETERIORATION

If turbidity has occurred, turbidity may be a sign of contamination.

WARNING AND PRECAUTIONS

1. The reagent is for *in vitro* diagnostic use. **Caution:** Do not pipette the solution by mouth. Avoid ingestion/contact.
2. Specimens should be considered infectious and handled appropriately.

STORAGE AND STABILITY

Store all reagent set at room temperature (15 - 30°C). The reagents are stable until the expiration date indicated on the label.

SPECIMEN COLLECTION

Freshly drawn serum is the specimen of choice and a 50 µl (0.05 ml) amount is required. Plasma from non-sodium containing anticoagulants (e.g., lithium, calcium, magnesium or heparin) is an acceptable alternative. Sodium is stable for at least 24 hours at room temperature and 2 weeks when refrigerated.^{1,2}

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1. Spectrophotometer.
2. Centrifuge.
3. Test tubes/rack.

PROCEDURE

Filtrate Preparation:

1. Label test tubes: blank, standard, control, patient, etc.
2. Pipette 1.0 ml of Filtrate Reagent to all tubes.
3. Add 50 µl of sample to all tubes and distilled water to the blank.
4. **Shake all tubes vigorously and mix continuously for 3 minutes.**
5. Centrifuge tubes at high speed (1,500G) for 10 minutes and test the supernatant fluids as described below, taking care not to disturb the protein precipitate.

COLOR DEVELOPMENT

1. Label test tubes corresponding to the above Filtrate tubes.
2. Pipette 1.0 ml of Acid Reagent to all tubes.
3. Add 50 µl of Supernatant to respective tubes and mix.
4. Add 50 µl of Color Reagent to all tubes and mix.
5. Zero spectrophotometer with distilled water at 550 nm.
6. Read and record absorbance of all tubes.

NOTE: The chemistry reaction of this procedure involves a reduction in absorbance, as opposed to the usual absorbance increase. The absorbance of blank should be higher than the test samples.

CALCULATIONS

Abs = Absorbance
S = Sample
STD = Standard

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$$\frac{(Abs. of Blank - Abs. of S)}{(Abs. of Blank - Abs. of STD)} \times Conc. of STD = Conc. of S$$

(mEq/L) (mEq/L)

SAMPLE CALCULATION

Assume the Standard with a sodium value of 150 mEq/L, gave an absorbance of 0.803 while the Sample and the Blank had absorbances of 0.880 and 1.406 respectively. The sodium concentration of the Sample may then be calculated as follows:

$$\frac{(1.406 - 0.880) \times 150}{(1.406 - 0.803)} = \frac{0.526 \times 150}{0.603} = 131 \text{ mEq/L}$$

PROCEDURAL LIMITATIONS

- When preparing filtrates, inadequate shaking or centrifugation will cause falsely lowered test results.
- Blood calcium, chloride, and potassium levels of up to 3 times normal reportedly exert no adverse influence on the procedure, phosphorus levels exceeding 5 times normal likewise present no problems.

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established sodium values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction, or procedural errors.

EXPECTED VALUES

135 - 155 mEq/L

PERFORMANCE CHARACTERISTICS

- Linearity:** 200 mEq/L.
- Sensitivity:** Based on an instrument resolution of $A = 0.001$, the present method has a sensitivity of 0.5 mEq/L.
- Comparison:** A comparison between this procedure and flame photometric analysis produced a regression equation of $Y = 0.69X + 4.5$ with a coefficient of correlation of 0.92.
- Precision Study**

Within Run		
Mean (mEq/L)	S.D.	C.V. %
146	7	5
127	4	3

Run to Run		
Mean (mEq/L)	S.D.	C.V. %
148	5	4
139	14	10

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High Technology, Inc

Type: Reagent
Method: Colorimetric
Form: Liquid
Linearity: up to 200 U/L
Sample: Urine

Kit Composition:
R1: 1 x 75 mL 0.2M MNP-glcNAc
(2-methoxy-4-(2-nitrovinyl)-phenyl-2-acetamido-2-deoxy-β-D-glucopyranoside)
R2: 1 x 15 mL 0.1M Citric Buffer
R3: 1 x 30 mL 0.5M Sodium Carbonate (Na₂CO₃) Buffer
Kit (540 Tests)

Intended use:
The N-acetyl-β-D-glucosaminidase (NAG) assay kit is for determination of NAG in urine samples.
FOR IN VITRO DIAGNOSTIC USE.

Storage and Handling:
The reagent is stable for 12 months from date of manufacture when stored at @ 2°-8°C.

Technical characteristics:
Precision:
The within-run and between-run precision has C.V. 5%.

INTERFERENCE:
Bilirubin: No interference up to 5 mg/dL
Hemoglobin: No interference
Albumin: No interference

Assay Procedure
An example of automated application:

R1+R2 150 μL
↓
□ ← Sample 10 μL
↓ 37°C, 5 min.
□ ← R3 (Sodium Carbonate Buffer) 50 μL
↓ 37°C, 20 Seconds.
Read at 505 nm wavelength
Note: Allow all reagents and specimens to warm to 37 degrees Celsius. Mix all reagents gently before using.

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INTENDED USE

For the quantitative determination of iron, total iron-binding capacity in human serum.

INTRODUCTION

The iron content of the human body may be divided into three classes: iron in storage, iron in use, and iron in transport. Iron in storage is reserved iron contained within the cells. Iron in use contained in hemoglobin, various enzymes, and several other types of proteins. Iron in transport is being moved to storage or is being removed from storage to be utilized in the formation of hemoglobin, etc. Iron in a free state is not only relatively insoluble, but it is toxic. Therefore, nearly all iron in the body is attached to some type of protein. It is of fundamental importance to note that a specimen should be analyzed for both iron and iron binding capacity because of the need for both values in the differential diagnosis of various types of anemia and liver diseases. For this reason, the current procedure is designed for the simultaneous determination of iron and iron binding capacity 1>2

Serum iron assays measure transport iron bound to the protein transferrin. Increase in serum iron levels may indicate increased erythrocyte destruction, decreased erythrocyte formation, increased absorption, or defects in storage capabilities. Decrease in serum iron levels may indicate iron deficiency or inability to retrieve storage iron. Iron binding capacity is usually increased in iron deficient anemia and decreased in hemochromatosis, malignancies, rheumatic fever, Hodgkin's diseases, and collagen vascular disease.^{1,2} Most successful iron methodologies remove iron from transferrin, reduce it to the ferrous state, bind it to a chromophore, and quantitate it by measuring the amount of color developed. The determination of iron binding capacity involves adding sufficient iron to saturate transferrin and then determining either the total amount of bound iron or the excess unbound iron. The latter procedure is applied to determine Iron Binding Capacity

PRINCIPLE

The iron in serum, is dissociated from its Fe (III) - transferrin complex by the addition of an acidic buffer containing hydroxylamine. This addition reduces the Fe (III) to Fe (II). The chromogenic agent, Ferene, forms a highly colored Fe (II) - complex that is measured photometrically at 560 nm. The unsaturated iron binding capacity (UIBC) is determined by adding Fe (II) iron to serum so that they bind to the unsaturated iron binding sites on transferrin. The excess Fe (II) ions are reacted with Ferrozine to form the color complex, which measured photometrically. The difference between the amount of Fe (II) added and the amount of Fe (II) measured represents the unsaturated iron binding. The total iron binding capacity (TIBC) is determined by adding the serum iron value to the UIBC value.

REAGENTS

1. Iron buffer reagent: Acetate buffer containing 220 mM Hydroxylamine hydrochloride, pH 4.5 with surfactant.
2. UIBC buffer reagent: Tris buffer 0.5 M, pH 8.0 with surfactant, sodium azide as preservative.
3. Iron color reagent: Ferrozine (16.6 mM) in Hydroxylamine hydrochloride.
4. Iron standard (500 ug/dl): 500 ug Ferrous chloride in Hydroxylamine hydrochloride.

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use.
 2. UIBC buffer contains sodium azide and may react with lead and copper plumbing to form highly explosive metal azides.
- On disposal, flush with a large volume of water to prevent

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High Technology, Inc

Total Iron binding capacity (TIBC) Reagent Set

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azide accumulation.

3. Avoid ingestion of reagent, as toxicity has not yet been determined.

4. Specimens should be considered infectious and handled appropriately.

STORAGE AND STABILITY

All the reagents and standard should be stored at room temperature (18-30°C).

REAGENT DETERIORATION

1. Appearance of turbidity, possible mold growth, or crystal formation that will not readily dissolve are signs of reagent deterioration.

2. Failure to obtain accurate results in the assay of control materials may indicate reagent deterioration.

SPECIMEN COLLECTION AND STORAGE

1. Fresh, unhemolyzed serum is the specimen of choice.

2. Serum should be separated as soon as clot has formed.

3. Heparinized plasma may be used but other anticoagulants should not be used to avoid possible iron contamination.

4. Serum iron is reported to be stable for four days at room temperature (18 - 30°C) and seven days at 2 - 8°C.

INTERFERING SUBSTANCES

1. Certain drugs and other substances are known to influence circulating iron levels. See Young, et al.⁴

2. Iron contained in hemoglobin does not react in this method; therefore, slight hemolysis will not interfere. However, gross hemolysis (pink or red specimens) will contribute to the absorbance measured at the wavelength used and should be avoided.¹

3. To make tubes, pipettes, etc. iron free, they must be washed with hot dilute 1:3 hydrochloric or nitric acid, followed by several rinsing with iron-free deionized or distilled water.

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1. Spectrophotometer capable of reading at 560 nm
2. Iron-free deionized water
3. Pipetting devices
4. Test tubes/rack
5. Timer
6. Heating bath/block

MANUAL PROCEDURE

Serum Iron

1. Label test tubes/cuvettes: "Blank", "Standard", "Control", "Sample", etc.

2. Add 2.5 ml Iron Buffer reagent to all tubes. 3. Add 0.5 ml (500µl) sample to respective tubes and mix. NOTE: Add 500 µl iron-free water to blank.

4. Zero spectrophotometer at 560 nm with the reagent blank.

5. Read and record the absorbances of all tubes (A₁ reading).

6. Add 0.05 ml (50 µl) iron color reagent to all the tubes. Mix.

7. Place all the tubes in the heating bath at 37°C for 10 minutes.

8. Zero the instrument at 560 nm with the reagent blank.

(Wavelength range: 520-560 nm)

9. Read and record absorbances of all the tubes (A₂ reading).

* HTI MULTI-PURPOSE CALIBRATOR MAYBE USED TO REPLACE STANDARD.

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$$\frac{0.07}{0.40}$$

$$0.175 \times 500 = 87.5 \text{ ug/dl}$$

UIBC (Unsaturated Iron-Binding Capacity)

1. Label test tubes/cuvettes, "Blank", "Standard", "Control", "Test", etc.
2. Add 2.0 ml UIBC buffer reagent to all tubes.
3. To "Blank" add 1.0 ml iron-free water. Mix.
4. To "Standard" add 0.5 ml (500 u/l) iron-free water plus 0.5 ml (500 u/l) standard. Mix.
5. To "Test" add 0.5 ml (500 u/l) respective sample plus 0.5 ml (500 u/l) iron standard. Mix.
6. Zero spectrophotometer at 560 nm with reagent blank.
7. Read and record the absorbance of all tubes (A₁ reading).
8. Add 0.05 ml (50 u/l) of Iron Color Reagent to all tubes. Mix.
9. Place all tubes in a heating bath at 37°C for ten (10) minutes.
10. Zero spectrophotometer at 560 nm with reagent blank.
11. Read and record the absorbance of all tubes. (A₂ reading).

HTI MULTI-PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

UIBC CALCULATIONS

$$\text{Conc. Stand.} - \frac{(A_{\text{2 Test}} - A_{\text{1 Test}})}{(A_{\text{2 Std}} - A_{\text{1 Std}})} \times \text{Conc. Stand.} = \text{UIBC (ug/dl)}$$

Example:

$$500 - \frac{(0.2 - 0.10)}{(0.4 - 0.00)} \times 500 = \text{UIBC (ug/dl)}$$

Therefore:

$$500 - (0.25 \times 500) = 375 \text{ ug/dl (UIBC)}$$

NOTE: The difference between A₁ Test and A₂ Test may sometimes be very small due to a high degree of unsaturation of transferrin with iron. The sample should be diluted with iron-free water and reassayed. Multiply the result by the dilution factor.

CALCULATIONS

TIBC (Total Iron-Binding Capacity):

Iron Level + UIBC = TIBC (ug/dl)

SI Unit Conversion: ug/dl x 0.179 = umol/L

CALIBRATION

The procedure is calibrated with iron standard (500 ug/dl) included in each kit.

QUALITY CONTROL

Serum controls with known normal and abnormal values should be run routinely to monitor the validity of the reaction.

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PERFORMANCE CHARACTERISTICS

1. Linearity, 500 ug/dl

Samples with values above 500 ug/dl must be diluted 1:1 with normal saline, re-assayed and result multiplied by two.

2. Sensitivity: Based on an instrument of A = 0.001, this procedure has a sensitivity of 2 ug/dl.

3. Comparison study: A study performed between this procedure and a similar Serum Iron procedure resulted in a coefficient of correlation of 0.98 with a regression equation of $y = 1.0x - 4.86$.

A study performed between this procedure and a similar UIBC procedure resulted in a coefficient of correlation of 0.97 with a regression equation of $y = 0.86x + 56$.

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The TBA assay is an enzyme cycling based assay utilizing the unique enzymatic property of 3- α hydroxy-steroid dehydrogenase (3- α HSD) that catalyzes the reversible reactions. The substrate bile acids is converted into 3-keto steroids in the forward reaction with concomitant conversion of Thio-NAD to Thio-NADH. In the reverse reaction, 3-keto steroids is converted back to bile acids with concomitant conversion of NADH to NAD. The amount of Thio-NADH formation in the cycling reaction is proportional to the amount of total bile acids in the sample. The formed Thio-NADH is detected at 405nm.

Reagent composition:

Reagent 1:
OVER 2.5 mmol / l
NBT 0.61 mmol / l
Diaphorase 625 U / L
Buffer pH 7.0 $\pm$ 0.1
Preservatives and stabilizers

Reagent 2:
Three- α -HSD 1250 U / L
Preservatives and stabilizers

Intended use

Total Bile Acids Assay Kit is intended for the in vitro quantitative determination of total bile acids (TBA) in human serum. Total bile acids are metabolized in the liver and serve as a marker for normal and abnormal liver function. Serum total bile acids are increased in patients with liver disease.

Product Packaging

Kit (400 Tests) R1: 2 x 60 mL
R2: 2 x 20 mL

Storage

The reagent is stable for 12 months from date of manufacture when stored 2°-8°C. Closed reagents can be used within 12 months from date of manufacture indicated on the packaging and labels.

Technical characteristics

Accuracy: the assay has excellent precision (%CV <5%) with no significant interference from lipemic or hemolytic samples.

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Type.: Reagent
Method: Immunoturbidimetric
Format: Liquid
Linearity: 10 - 2,000 IU/mL (or value of highest calibration point)
Sample: Serum, Plasma

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KIT COMPOSITION

R1: 2 X 20 mL buffer
R2: 1 X 6.6 mL antibody
Tests / Kit: 220 (Hitachi 917)
(Can Vary Based On Analyzer Model)

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INTENDED USE

The K-ASSAY® Total IgE is an in vitro diagnostic reagent for the quantitative determination of circulating total IgE in human serum or plasma by immunoturbidimetric assay for use as an aid in the clinical diagnosis of IgE mediated allergic disorders in conjunction with other clinical findings.
FOR IN VITRO DIAGNOSTIC USE.

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PERFORMANCE

Precision

The precision for the K-ASSAY® Total IgE assay was determined using packaged reagents, pooled human serum, and a Hitachi 917 analyzer.

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Accuracy / Correlation

A comparison of the K-ASSAY® Total IgE assay and Company X's Total IgE reagent was performed with the following results:

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Assay Range

10 - 2,000 IU/mL (or value of highest calibration point)

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INTERFERENCE

Bilirubin F :	No interference up to 18.5 mg/dL
Bilirubin C :	No interference up to 21.6 mg/dL
Hemoglobin :	No interference up to 467 mg/dL

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RF : No interference up to 500 IU/mL
Sodium Citrate : No interference up to 15 mg/mL
Sodium Heparin : No interference up to 150 IU/mL
EDTA-2Na : No interference up to 1.5 mg/mL
Chyle (Formazine Turbidity) : No interference up to 2,120

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Type: Reagent
Method: Enzyme
Linearity: up to 1000 $\mu\text{mol/L}$
Sample: serum or plasma

Kit Composition:
Sample: 4 μL
R1: 2x27ml 0.3M Thio-NAD in Buffer
(3-(N-morpholino)propanesulfonic acid)
R2: 2x9ml 0.1M Acetoacetyl-coA thiolase (Enzyme)

Intended Use:
The Total Ketone Bodies is an in vitro assay for the quantitative determination of total ketone bodies [acetoacetone (AcAc) + 3-hydroxybutyrate (3-HB)] in serum or plasma.
FOR IN VITRO DIAGNOSTIC USE.

Storage and handling:
All reagents must be stored at 2-10°C.

Assay Procedure
An example of automated application:
Sample 4 μL
↓
□ ← R1 (Thio-NAD) 270 μL
↓ 37°C, 5 min.
□ ← R2 (Enzyme) 90 μL
↓ 37°C, 5 min.
Rate, 400 nM (main) / 600 nM (sub)
Note: Mix all reagents gently before using.

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TYPE: Reagent

METHOD: Immunoturbidimetric

FORMAT: Liquid (Lyophilized Calibrator)

LINEARITY: 3 - 14 mg/dL

SAMPLE: Plasma

KIT COMPOSITION

R1: 4 X 20 mL buffer

R2: 2 X 10 mL antibody

Calibrator: 1 x 1 mL (lyophilized)

Tests / Kit: 250 (Hitachi 717)

(Can Vary Based On Analyzer Model)

INTENDED USE:

For the quantitative determination of human plasminogen in human plasma by immunoturbidimetric assay.

FOR IN VITRO DIAGNOSTIC USE

REAGENT PREPARATION:

Reagent R1 and R2 are ready to use and do not require reconstitution.

The calibrator is lyophilized and requires reconstitution according to the package insert instructions.

STORAGE AND HANDLING:

All reagents should be stored refrigerated (2-8°C). Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for one year from the date of manufacture, as indicated by the expiration date on the package and bottle labels.

Opened R1 and R2 can be used for one month if stored at 2-8°C.

Reconstituted plasminogen calibrator can be used for 1 week after reconstituting if stored at 2-8°C.

PERFORMANCE:

Precision

The precision for Plasminogen assay was determined using packaged reagents, pooled human plasma, and a Hitachi 717 analyzer.

Precision - Within Run

Sample I

N = 20

Mean = 5.33

S.D. = 0.086

C.V. = 1.62%

Precision - Between Runs

Sample I

N = 10

Mean = 5.31

S.D. = 0.074

C.V. = 1.39%

Sample II

N = 20

Mean = 11.4

S.D. = 0.122

C.V. = 1.07%

Sample II

N = 10

Mean = 11.5

S.D. = 0.099

C.V. = 0.86%

Accuracy / Correlation

A comparison of Plasminogen and a Binding Site Plasminogen RID Kit was performed using a Hitachi 717. The test results provided the following data: $y = 0.995x + 0.309$

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INTENDED USE

For the qualitative and semi-quantitative measurement of RF in human serum.

INTRODUCTION

Rheumatoid factors (RF) are antibodies directed against antigenic sites in the Fc fragment of human and animal IgG. Their frequent occurrence in rheumatoid arthritis makes them useful for diagnosis and monitoring of the disease.^{1,2} One method used for rheumatoid factor detection is based on the ability of rheumatoid arthritis sera to agglutinate sensitized sheep red cells, as observed by Waaler³ and Rose.⁴ A more sensitive reagent consisting of biologically inert latex beads coated with human gamma globulin was later described by Singer and Plotz.⁵ The RF kit is based on the principle of latex agglutination assay by Singer and Plotz.⁵ The major advantage of this method is rapid performance (3 minute reaction time) and lack of heterophile antibody interference.

PRINCIPLE

The RF reagent is a suspension of polystyrene latex particles sensitized with specially prepared human IgG. The reagent is based on an immunological reaction between human IgG bound to biologically inert latex particles and rheumatoid factors in the test specimen. When serum containing rheumatoid factors is mixed with the latex reagent, visible agglutination occurs. The RF latex reagent sensitivity has been adjusted to detect a minimum of 8 IU/mL of rheumatoid factors according with the WHO International Standard without previous sample dilution.

REAGENTS

1. **RF Latex Reagent:** A suspension of uniform polystyrene particles coated with IgG (human) in glycine buffer, pH 8.2; reagent sensitivity is standardized with the World Health Organization RF Standard, MIX WELL BEFORE USING.
2. **RF Positive Control Serum:** A stabilized, prediluted human serum containing at least 8 IU/mL of RF.
3. **RF Negative Control Serum:** A stabilized, prediluted human serum containing less than 8 IU/mL of RF.
4. **Glycine-Saline Buffer (20x):** pH 8.2 ± 0.1M glycine and 0.15M NaCl.

NOTE

1. Dilute buffer following instructions on the label before using.
2. All reagents contain 0.1% (w/v) sodium azide as a preservative.
3. The human sera used in the controls has been tested and found negative for HbsAg and HIV, however, careful handling is always recommended.
4. Store all reagents at 2 - 8°C. DO NOT FREEZE.

WARNINGS AND PRECAUTIONS

1. Reagents containing sodium azide may combine with copper and lead plumbing to form highly explosive metal azides. Dispose of reagents by flushing with large amounts of water to prevent azide buildup.

For In Vitro diagnostic use. Positive and negative controls prepared using human sera were found negative for hepatitis B surface antigen (HBsAg) as required by FDA; however, handle controls as if potentially infectious.

REAGENT STORAGE AND STABILITY

1. Reagents are stable until stated expiration date on bottle label when stored refrigerated (2 - 8°C).
2. DO NOT FREEZE.
3. The RF latex reagent, once shaken must be uniform without visible clumping. When stored refrigerated, a slight sedimentation may occur and should be considered normal.
4. Do not use the latex reagent or controls if they become contaminated.

SPECIMEN COLLECTION AND STORAGE

1. Use fresh serum collected by centrifuging clotted blood.
2. If the test cannot be carried out on the same day, the serum may be stored between 2 - 8°C for no longer than 72 hours after collection.
3. For longer periods the sample must be frozen.
4. As in all serological tests, hemolytic or contaminated serum must not be used.
5. Do not use plasma.

MATERIALS AND REAGENTS PROVIDED

1. RF Latex Reagent
2. RF Positive Control
3. RF Negative Control
4. Glycine-Saline Buffer
5. Reaction Slide
6. Pipette/Stir Sticks

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Qualitative Test

1. Bring reagents and specimens to room temperature before use.
2. Place one drop (50 µl) of the RF Positive Control on field #1 of the reaction slide. Place one drop (50 µl) of the RF Negative Control on field #2. The remaining fields are used for test specimens. Using pipettes provided, place one drop of the undiluted specimens on successive fields. Retain Pipette/Stir Sticks for mixing step.
3. Gently resuspend the RF Latex Reagent and add one drop to each test field. Use pipette/Stir Stick to spread reaction mixture over entire test field.
4. Rotate the slide manually or with a mechanical rotor at 80-100 rpm for 2 minutes and read immediately under direct light.
5. Presence of agglutination of the latex particle is a positive result (see figure 1). Agglutination indicates a RF concentration of equal or more than 8 IU/ml. Sera with positive agglutination should be run again with the Quantitative Test.

Quantitative Test

1. Bring reagents and specimens to room temperature before use.
2. Using the Glycine-Saline Buffer, dilute the specimens 1:2, 1:4, 1:8, 1:16, 1:32 or as needed.
3. Place one drop (50 µl) each of negative and positive controls on two slide rings. Place on drop (50 µl) of each dilution on successive fields of the reaction slide.
4. Gently resuspend the RF Latex Reagent and add one drop to the reaction slide.
5. Gently resuspend the RF Latex Reagent and add one drop to each test field. Use Pipette/Stir Stick to spread reaction mixture over entire field.
6. Rotate the slide for 2 minutes and read immediately under direct light.

QUALITY CONTROL

1. RF Positive and Negative Control should be included in each test batch.
2. Acceptable performance is indicated when a uniform milky suspension with no agglutination is observed with the RF Negative Control and agglutination with large aggregates is observed with the RF Positive Control (Figure 1).

INTERPRETATION

Negative Result

A negative reaction is indicated by a uniform milky suspension with no agglutination observed with the RF Negative Control.

Positive Result: A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared to the RF Negative and Positive Controls (Fig. 1).

Semi-quantitative Test

The titer of the serum is the reciprocal of the highest dilution, which exhibits a positive reaction.

An estimate of the RF concentration in the specimen can be expressed in IU/ml by using the following equation:

$$\text{IU/ml of specimen} = \text{IU/ml control} \times \text{specimen titer}$$

LIMITATIONS

1. Results should be read two (2) minutes after the mixing or the reagent on the slide.
2. Existence of prozone at high titers has not been encountered.
3. Increased levels of RF may be found in some diseases other than rheumatoid arthritis such as infectious mononucleosis, sarcoidosis, lupus erythematosus, Sjogren's syndrome, etc.
4. Certain patients with rheumatoid arthritis will not have the RF present in their serum.

EXPECTED VALUES

1. The diagnosis of rheumatoid arthritis is based largely on clinical examination, but laboratory tests are useful to support the clinical diagnosis and to evaluate the severity and course of the disease in the individual patient. One

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of the most useful clinical markers for rheumatoid arthritis is rheumatoid factor in serum. Rheumatoid factor is a term used to describe a variety of antibodies or immune complexes or both, that occur with rheumatoid arthritis as well as in a variety of other diseases.⁸

2. Different studies have shown positive serological reactions for rheumatoid factor in as high as 90% of patients with rheumatoid arthritis compared with less than 5% in control groups.¹

PERFORMANCE

1. Sensitivity: 8 IU/ml or above.

2. Comparison:

A. Qualitative Results: The RF Latex Reagent was evaluated on a total of 75 samples from the hospital patients. The qualitative test was evaluated by comparison with a commercially available latex agglutination test. This study demonstrated a 96% agreement between these tests. The discrepancy results were obtained in samples with titers near the limit of sensitivity of the reagents.

B. Semi-quantitative Results: A panel of 10 known RF positive serum samples was quantitated on three consecutive days. The results of the study indicated that RF Latex Reagent has 100% precision.

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High Technology, Inc.

INTENDED USE

C-Reactive Protein (CRP) latex slide test is used for the qualitative and semi-quantitative measurement of C-reactive protein (CRP) in human serum.

INTRODUCTION

C-Reactive Protein (CRP), the classic acute-phase of human serum, is synthesized by hepatocytes. Normally, it is present only in trace amounts in serum, but it can increase by as much as 1000-fold in response to injury or infection. The clinical measurement of CRP in serum, therefore, appears to be a valuable screening test for organic disease and a sensitive index of disease activity in inflammatory, infective, and ischemic conditions.^{1,2} MacLeod and Avery found that antibody produced against purified CRP provided a more sensitive test than the C-polysaccharide assay.³ Since that time a number of immunological assays have been devised to measure CRP such as capillary precipitation, double immunodiffusion, and radical immunodiffusion.^{4,5} The CRP reagent kit is based on the principle of the latex agglutination assay described by Singer and Plotz.⁶ The major advantage of this method is the rapid three (3) minute reaction time.

PRINCIPLE

The CRP reagent kit is based on an immunological reaction between CRP antisera bound to biologically inert latex particles and CRP in the test specimen. When serum containing greater than 0.8 mg/dl CRP is mixed with the latex reagent, visible agglutination occurs.

REAGENTS

1. CRP Latex Reagent: A suspension of uniform polystyrene particles coated with monospecific antihuman CRP (goat) in glycine buffer, pH 8.8 ± 0.5; reagent sensitivity adjusted to approximately 0.8 mg/dl. MIX WELL BEFORE USING.
2. CRP Positive Control Serum: A stabilized prediluted human serum containing more than 0.8 mg/dl CRP.
3. CRP Negative Control Serum: A stabilized prediluted human serum non-reactive with the test reagent.
4. Glycine-Saline Buffer (20x): pH 8.2 ± 0.1. A diluent containing 0.1M glycine and 0.15M NaCl. Dilute glycine-saline buffer before using by adding 1 part glycine-saline buffer solution to 19 parts distilled water. All reagents contain 0.1% (w/v) sodium azide as preservative.

WARNINGS AND PRECAUTIONS

1. Reagents containing sodium azide may combine with copper and lead plumbing to form highly explosive metal azides. Dispose of reagents by flushing with large amounts of water to prevent azide buildup.
2. For in vitro diagnostic use.
3. Positive and negative controls prepared using human sera found negative for hepatitis B surface antigen (HBsAg); however, handle controls as if potentially infectious.

REAGENT STORAGE AND STABILITY

1. Reagents are stable until stated expiration date on bottle label when stored refrigerated (2 - 8°C).
2. DO NOT FREEZE.
3. The CRP latex reagent, once shaken must be uniform without visible clumping. When stored refrigerated, slight sedimentation may occur and should be considered normal.
4. Do not use the latex reagent or controls if they become contaminated.

SPECIMEN COLLECTION AND STORAGE

1. Use fresh serum collected by centrifuging clotted blood.
2. If the test cannot be carried out on the same day, the serum

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High Technology, Inc

may be stored between 2 - 8°C for no longer than 72 hours after collection.

- For longer periods, the sample must be frozen.
- As in all serological tests, hemolytic or contaminated serum must not be used.
- Do not use plasma.

MATERIALS AND REAGENTS PROVIDED

- CRP Latex Reagent
- CRP Positive Control
- CRP Negative Control
- Glycine-Saline Buffer
- Reaction Slide
- Pipette/Stir Sticks

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PROCEDURE

QUALITATIVE TEST:

- Bring reagents and specimens to room temperature before use.
- Place one drop (50 µl) of the CRP Positive Control on field #1 of the reaction slide. Place one drop (50 µl) of the CRP Negative Control on field #2. The remaining fields are used for test specimens. Using pipettes provided, place one drop of the undiluted specimens on successive fields. Retain pipette/stir sticks for mixing step.
- Gently resuspend the CRP Latex Reagent and add one drop to each test field. Use pipette/stir stick to spread reaction mixture over the entire test field.
- Rotate the slide for three (3) minutes and read immediately under direct light.

SEMI-QUANTITATIVE TEST:

- Set up at least five test tubes: 1:2, 1:4, 1:8, 1:16, etc.
- Dilute sample according to dilution factor on each test tube with diluted saline solution.
- Place one drop of each of positive and negative controls onto separate slide fields. Place one drop of each dilution on successive fields of the reaction slides.
- Gently resuspend the CRP Latex Reagent and add one drop to each test field.
- Mix well with the flat end of the pipette. Gently rock the slide for three (3) minutes and read immediately under direct light. The titer of the serum is the reciprocal of the highest dilution exhibiting a positive reaction multiplied by the concentration of the positive control.

concentration reciprocal of last
mg/dL of serum = of positive × dilution showing
control a positive result

Dilution	Reciprocal	Conc. of serum (mg/dL)
1:2	2	1.6
1:4	4	3.2
1:8	8	6.4
etc.		

QUALITY CONTROL

- CRP Positive and Negative Control should be included in each test batch.
- Acceptable performance is indicated when a uniform milky suspension with no agglutination is observed with the CRP Negative Control and agglutination with large aggregates is observed with the CRP Positive Control.

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High Technology, Inc.

Negative Result: A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the CRP Negative Control.

Positive Result: A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared to the CRP Negative Control (Figure 1).


Figure 1.

EXPECTED VALUES

1. CRP in healthy individuals is approximately 0.02-1.35mg/dl. The mean value in adults is 0.047mg/dl.
2. A weak positive correlation was found between CRP and age.
3. It is important to determine the level of CRP for monitoring patient progress. This is due to (1) the concentration of CRP is an index of tissue damage incurred and (2) increasing or decreasing levels of CRP (e.g. daily) indicate the progress of inflammatory process.

LIMITATIONS

1. Reaction time is critical. If reaction time exceeds three (3) minutes, drying of the reaction mixture may cause false positive results.
2. Freezing the CRP Latex Reagent will result in spontaneous agglutination.
3. Intensity of agglutination is not necessarily indicative of relative CRP concentration; therefore, screening reactions should not be graded.
4. A false negative can be attributed to a prozone phenomena (antigen excess). It is recommended, therefore, to check all negative sera by retesting at a 1:10 dilution with glycine buffer.

PERFORMANCE

1. **Sensitivity:** greater than 0.8 mg/dl
2. **Linearity:** 25 mg/dl.

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Kit Composition:
R1: 1 X 30 mL buffer
R2: 1 X 20 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model And Calibrator Set Used)

Intended Use:
For the quantitative determination of C-reactive protein in serum by latex particle enhanced immunoturbidimetric assay. Measurement of C-reactive protein aids in evaluation of the amount of injury to body tissues.
FOR IN VITRO DIAGNOSTIC USE.

Reagent Preparation:
The reagents are ready to use and do not require reconstitution. Mix gently before using.

Storage and Handling:
All reagents should be at 2-8°C and protected from light. Unopened reagents can be used for 18 months from the date of manufacture as indicated on the expiration date on the package and bottle labels. Once the reagent vial has been opened, store tightly capped at 2-8°C and use within 1 month. It is recommended that R-2 be gently inverted once a week after it is opened.

Performance:
Precision
Samples tested were commercial human CRP control serum

Precision - Within Run		
Sample I	Sample II	Sample III
N=20	N=20	N=20
Mean = 2.74 mg/L	Mean = 0.274 mg/dL	Mean = 8.64 mg/L
	or	
Mean = 0.864 mg/dL	Mean = 17.43 mg/L	Mean = 1.743 mg/dL
S.D. = 0.034	S.D. = 0.068	S.D. = 0.135
C.V. = 1.24%	C.V. = 0.78%	C.V. = 0.77%

Precision - Between Runs		
Sample I	Sample II	Sample III
N=10	N=10	N=10
Mean = 2.65 mg/L	Mean = 8.55 mg/L	Mean = 17.13 mg/L
	or	
Mean = 0.265 mg/dL	Mean = 0.855 mg/dL	Mean = 1.713 mg/dL
S.D. = 0.111	S.D. = 0.092	S.D. = 0.204
C.V. = 4.2%	C.V. = 1.08%	C.V. = 1.19%

Accuracy / Correlation

$y = 1.034x - 0.173$
 $r = 0.999$
 x = Company A's latex CRP nephelometric assay
 y = CRP (Z)

Assay Range
With CRP (2) Calibrator Set A
0.1 - 20 mg/L
or
0.01 - 2.0 mg/dL

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High Technology, Inc

For the In Vitro Quantitative Determination of Salicylate in Serum

Summary and Principle

Salicylate is a common drug used in many formulations due to its analgesic and anti-inflammatory properties. Salicylate intoxication occurs often in children due to its accessibility, in chronic patients with regular need of medication or in patients who are taking combination prescription and non-prescription dosages. Overdosage of salicylate is also associated with suicide attempts in adolescents and adults. 1 Serum concentrations exceeding 600 mg/L are usually lethal. 2,3,4 Salicylate overdose can cause metabolic acidosis with a high anion gap, gastrointestinal and central nervous system disturbances, as well as encephalopathy and renal failure. 5 Immediate detection and treatment of this medical emergency is essential. Therefore, a method for the rapid and accurate determination of salicylate is needed. Enzymatic methods used for the quantitation of salicylate, as opposed to traditional chemical tests, 6 provide specificity and simplicity. In addition, the Stanbio Salicylate Reagent offers a rapid and convenient method for salicylate quantitation.

The present enzymatic assay for the determination of salicylate depends upon the conversion of salicylate in the presence of NADH by Salicylate hydroxylase (SH) to catechol and NAD.

Salicylate
hydroxylase

The conversion of NADH to NAD is measured by the decrease in absorbance at 340 nm. The decrease is proportional to the concentration of salicylate present in the sample.

Reagents

Reagent (R1), Ref. No. 2411 Contains 0.3 mM NADH

Reagent (R2), Ref. No. 2412 Contains 20,000 U/L Salicylate Hydroxylase

Standard - 300 mg/L, Ref. No. 2413 Contains Salicylate in an aqueous buffer

Precautions

For in vitro diagnostic use only. No special precautions are needed with these reagents. However, general care in reagent handling is recommended. Handle serum and blood specimens as potentially infectious samples and follow the guidelines established by the Centers for Disease Control (CDC), Atlanta, Georgia, for blood collection and handling.

Reagent Preparation:

Reagents are ready to use.

Reagent Storage and Stability:

The reagents are stable for 24 months when stored at 2-8 °C when kept separated. The aqueous Salicylate Standard is stable at 2-8 °C for 24 months after manufacture. Do Not Freeze.

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Serum is the recommended specimen.

Sample Stability: Serum containing salicylate can be stored at 2-8 °C for one week.

Interfering Substances: The following substances, when added at a concentration of 500 mg/L to a serum containing 280 mg/L of salicylate, showed no interference.

Substance % Recovery	
Acetaminophen	103
Acetylsalicylic acid	105
Bilirubin	103
EDTA disodium	102
Ibuprofen	100
α-Ketobutyric acid	100
Methyl salicylate	100
Phenol	103
Salicylamide	102
Sodium benzoate	104
Sodium oxalate	100
Theophylline	99
Uric acid	104

Ethanol was tested at a concentration of 5% and found not to interfere.

Heparin at a concentration of 14.3 u/L and hemoglobin at OD540 of 2.0 were tested and found not to interfere. Additionally, when sodium citrate, sodium fluoride and oxalic acid at concentrations of 3.5 mg/mL, 2.5 mg/mL, and 2.0 mg/mL, respectively, were added to the same serum pool containing 280 mg/L of salicylate, no interference was observed.

However, in serum containing 280 mg/L salicylate, 500 mg/L concentrations of the drugs p-aminosalicylic acid, 2,5-dihydroxybenzoic acid and 2,5-dihydroxyphenylacetate gave biased values of +115%, +23% and +19% respectively.

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High Technology, Inc

The Salicylate Liqui-UV® Test can be used in either 5 minute procedures below. Any spectrophotometer that reads absorbance to the third decimal place may be used.

Conditions: Wavelength	340 nm
Temperature	37°C
Mode	End Point
Sample Volume	50 µL
Reagent Volume	2.04 mL
Sample to reagent ratio	1/42

Procedure #1

- 1) Into four 3 mL cuvettes add 2.0 mL of Reagent A at 37 °C.
- 2) Add 50 µL of serum to cuvettes 1 and 2.
- 3) Add 50 µL of Salicylate Standard to cuvettes 3 and 4.
- 4) Add 40 µL of water to cuvette 1 (serum blank) and cuvette 3 (std. blank).
- 5) Add 40 µL of Reagent B to cuvettes 2 and 4.
- 6) Measure the absorbance of all cuvettes 5 minutes after the addition of Reagent B.
- 7) Subtract the absorbance of cuvette 2 (ODserum) from the absorbance of cuvette 1 (ODserum blank).
- 8) Subtract the absorbance of cuvette 4 (ODstd) from the absorbance of cuvette 3 (ODstd blank).

Calculations

ODserum blk - ODserum x 300 mg/L x dil of serum = Salicylate mg/L

ODstd blk - ODstd

Procedure #2

- 1) Into two 3 mL cuvettes add 2.0 mL of Reagent A at 37 °C.
- 2) To cuvette 1 add 50µL of serum and measure the initial OD340 (TO).
- 3) To cuvette 2, add 50 µL of Salicylate Standard and measure the initial OD340 (Tstd).
- 4) To cuvette 1, add 40 µL of Reagent B and measure the final OD340 (TF) at 5 minutes.
- 5) To cuvette 2, add 40 µL of Reagent B and measure the final OD340 (TF std) at 5 minutes.
- 6) Calculate the ΔODS. (OD To x 0.981) - OD TF

Calculations

ΔOD 5 serum x 300 mg/L x dil of serum = Salicylate mg/L

Quality Control

In order to assure consistent performance, it is recommended that a normal and an abnormal serum control be assayed with each run.

Expected Values

The quantitation of salicylate is important in cases of drug overdose. Serum concentrations above 300 mg/L are generally toxic and concentrations exceeding 600 mg/L are usually lethal.

Performance Characteristics

Linearity: The procedure described is linear to 1000 mg/L salicylate. For higher concentrations, dilute the samples with distilled water. Repeat the assay and multiply the results by the dilution factor.

Sensitivity: Concentrations of salicylate of about 60, 90 and 100 mg/L can be clearly distinguished at the 99.9% confidence limit with this method.

Accuracy: Analytical recovery of salicylate added to human serum pools gave a mean recovery of 96% from the following individual values.

% Recovery

Pool 1 (100 mg/L)	99
Pool 2 (200 mg/L)	95
Pool 3 (300 mg/L)	97
Pool 4 (400 mg/L)	94

Comparison of 30 serum samples between the Stanbio Salicylate Liqui-UV® and the DCL Enzymatic Salicylate Assay Reagent Kit (Diagnostic Chemicals Limited) showed a correlation coefficient of r=0.9939. The comparison of the two methods produced the following correlation equation: Stanbio=0.914DCL + 4.047 mg/L, Std. Error: Syx=13.8 mg/L.

Precision: Studies were conducted using 3 serum pools containing 50, 150 and 300 mg/L of salicylate. The following results are an average of 18

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KIT COMPOSITION:

Option 1

R1: 3 X 20 mL buffer
R2: 2 X 10 mL antibody
Tests / Kit: 200 (Hitachi 717)
(Can Vary Based On Analyzer Model)

Option 2

Note: Large size kit, available by special order.
R1: 3 X 60 mL buffer
R2: 1 X 60 mL antibody

Tests / Kit: 600 (Hitachi 717)
(Can Vary Based On Analyzer Model)

INTENDED USE:

Transferrin is intended for the quantitative determination of human transferrin in human serum by immunoturbidimetric assay.
FOR IN VITRO DIAGNOSTIC USE.

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PERFORMANCE:

Precision:

The precision for the Transferrin was determined using packaged reagents, pooled human serum, and a Hitachi 717 analyzer.

Precision - Within Run

Sample I	Sample II	Sample III
N = 20	N = 20	N = 20
Mean = 70.7	Mean = 300.75	Mean = 526.4
S.D. = 0.92 S.D. = 5.02	S.D. = 5.77	
C.V. = 1.3%	C.V. = 1.67%	C.V. = 1.097%

Precision - Between Runs

Sample I	Sample II
N = 8	N = 8
Mean = 175.38	Mean = 340.9
S.D. = 1.6	S.D. = 1.46
C.V. = 0.91%	C.V. = 0.4%

Accuracy / Correlation

A comparison of the Transferrin to an Incstar Transferrin Test Kit was performed using a Hitachi 717. The test results provided the following data:

$y = 1.064x - 50.124$
 $r = 0.982$
 $n = 60$
 $x = \text{Incstar Transferrin Test Kit}$
 $y = \text{Transferrin}$
 $x_{\min} = 148$ $y_{\min} = 110$
 $x_{\max} = 440$ $y_{\max} = 426$
 $x_{\text{mean}} = 246$ $y_{\text{mean}} = 212$

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INTENDED USE

For the *In Vitro* quantitative determination of Triglycerides in serum or plasma.

INTRODUCTION

Triglycerides are esters of fatty acids and are hydrolyzed to glycerol and free fatty acids. Triglyceride determinations when performed in conjunction with other lipid assays are useful in the diagnosis of primary and secondary hyperlipoproteinemia. They are also of interest in following the course of diabetes mellitus, nephrosis, biliary obstruction and various metabolic abnormalities due to endocrine disturbances.

Standard methods for the measurement of triglyceride concentrations have involved either enzymatic or alkaline hydrolysis to liberate glycerol. This formulation makes use of the enzymatic hydrolysis and quantification since it is specific and not subject to interference by phospholipids.

PRINCIPLE

The enzymatic reaction sequence employed in the assay of Triglycerides is as follows:

Triglycerides + Lipase $\rightarrow$ Glycerol + fatty Acids

Glycerol + ATP Glycerol-Kinase $\rightarrow$ Glycerol-1-phosphate + ADP

Glycerol-1-Phosphate + O_2 -NAD $\rightarrow$ DAP + H_2O_2

H_2O_2 + 4-AAA + DHBS Peroxidase $\rightarrow$ Quinoneimine Dye + $2H_2O$

The present procedure involves hydrolysis of triglycerides by lipase. The glycerol concentration is then determined by enzymatic assay coupled with Trinder reaction that terminates in the formation of a quinoneimine dye. The amount of the dye formed, determined by its absorption at 520 nm, is directly proportional to the concentration of triglycerides in the samples.^{1,2}

REAGENT COMPOSITION

1. Triglyceride Liquid reagent contains the following:

ATP 0.5 mmol/L, Magnesium acetate 12 mmol/L, 4-Chlorophenol 3.5 mmol/L, 4-Aminophenazon 0.3 mmol/L, Glycerol Phosphate Oxidase >4500 U/L, Lipase $>200,000$ U/L, Glycerol Kinase >250 U/L, Peroxidase $>2,000$ U/L, Buffer (pH 7.4) 50 mmol/L, surfactants, stabilizers, and preservatives.

2. Triglyceride standard contains glycerol with surfactant to yield 200 mg/dl triglycerides as triolein. Sodium azide 0.1% is added as a preservative.

WARNINGS AND PRECAUTIONS

1. For *In Vitro* diagnostic use.

2. Avoid ingestion of reagent as toxicity has not yet been determined.

3. All specimens and controls should be considered infectious and handled appropriately.

4. Reagent and standard contain sodium azide as a preservative. This may react with copper or lead plumbing to form explosive metal azides. Upon disposal, flush with large amount of water to prevent azide build up.

REAGENT PREPARATION

Triglyceride reagent and standard are provided in a ready-to-use form.

No preparation is necessary.

STORAGE AND STABILITY

Both the Triglyceride reagent and standard must be stored at $2 - 8^\circ C$. The reagent may be used until the expiration date indicated on the package label when stored as directed. Protect from direct light. Avoid microbial contamination.

REAGENT DETERIORATION

The reagent should be discarded if:

1. The initial absorbance of the reagent against water is greater than 0.500 when measured at 520 nm.

2. The reagent fails to meet linearity claims or fails to recover stated values.

Note: A yellow or pink coloration is normal.

3. The reagent is turbid or displays evidence of bacterial contamination.

SPECIMEN COLLECTION

1. Fresh, clear, non-hemolyzed serum from fasting patients is recommended.

2. Triglycerides in serum appear stable for three (3) days when stored at $2 - 8^\circ C$.

3. Prolonged storage of the samples at room temperature is not recommended since other glycerol containing compounds may hydrolyze, releasing free glycerol.

4. Blood collection devices lubricated with glycerin (glycerol) should not be used.

INTERFERENCES

Glycerol in rubber stoppers or in contaminated glassware will elevate triglyceride levels. Lipemic or grossly icteric samples will cause falsely elevated results consequently a patient blank should be run. Samples with gross hemolysis or high bilirubin values will produce falsely elevated triglyceride values. A number of drugs and substances affect the measurement of triglyceride.³

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High Technology, Inc

1. Spectrophotometer capable of measuring absorbances at 500-550 nm
2. Test tubes and rack
3. Accurate pipetting devices
4. Constant temperature incubator set at 37° C
5. Timer

GENERAL INSTRUCTIONS

The reagent for triglycerides is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

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Label tubes: "blank", "standard", "control", "patient", etc.

Pipette 1.0 ml of reagent into all tubes.

Place all tubes in a 37° C heating block for at least 4 minutes.

Add 0.010 ml (10 µl) of sample to respective tubes and mix.

Incubate all tubes for five (5) minutes at 37° C.

Zero spectrophotometer at 520 nm with reagent blank

(Wavelength range 500-550)

Read and record absorbances of all tubes.

MANUAL PROCEDURE

Label tubes: "blank", "standard", "control", "patient", etc.

Pipette 1.0 ml of reagent into all tubes.

Place all tubes in a 37° C heating block for at least 4 minutes.

Add 0.010 ml (10 µl) of sample to respective tubes and mix.

Incubate all tubes for five (5) minutes at 37° C.

Zero spectrophotometer at 520 nm with reagent blank

(Wavelength range 500-550)

Read and record absorbances of all tubes.

Note: Final color is stable for sixty (60) minutes at room temperature.

* HTI - MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

PROCEDURAL LIMITATIONS

The reagent is linear to 1000 mg/dl (11.4 mmol/L), specimens above this limit must be diluted 1:1 with water, reassayed and multiplied the results by two to compensate for the dilution.

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established triglyceride values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedural errors.

CALCULATIONS

Triglycerides results are expressed as mg/dl

A = Absorbance

$$\frac{A(\text{patient})}{A(\text{standard})} \times \text{Concentration of patient} = \text{Concentration of patient (mg/dl)}$$

Example:

$$\frac{0.24}{0.31} \times 200 = 154.8 \text{ mg/dl}$$

NOTE:

To obtain the results in SI units (mmol/L) multiply the result in mg/dl by 0.0113.

EXPECTED VALUES

36 - 165 mg/dl

It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE CHARACTERISTICS Linearity 1000 mg/dl

Sensitivity: Based on an instrument resolution of A = 0.001, this procedure has a sensitivity of 1.3 mg/dl.

Comparison: A group of 91 sera ranging in Triglyceride values from 12 to 1030 mg/dl were assayed by this method and by a similar commercially available reagent. Comparison of the results yielded a correlation coefficient of 0.997 and the regression equation was $y = 0.946x + 5.373$. (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T.)

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INTENDED USE


For manual or automated quantitative determination of Carbon Dioxide content in human plasma or serum at 340 nm.

INTRODUCTION

Carbon Dioxide (CO₂) in serum or plasma exists primarily as dissolved CO₂ and bicarbonate anion (HCO₃⁻). The plasma CO₂ content is decreased in metabolic acidosis and respiratory alkalosis, whereas the level is increased in metabolic alkalosis and respiratory acidosis. In pathologic conditions such as in diabetes mellitus, glomerulonephritis, pyloric obstruction, diarrhea, etc., acidosis or alkalosis could be anticipated. Therefore, determination of plasma CO₂ content as part of an electrolyte profile can help establish to a degree, the anticipated change in the patient.

Plasma CO₂ content can be measured manometrically or using an ion selective electrode or by spectrophotometric procedures. This CO₂ reagent measures CO₂ content enzymatically and the procedure is a modification of the method of Forrester et al.

PRINCIPLES OF THE PROCEDURE

Phosphoenolpyruvate carboxylase (PEPC) catalyzes the reaction between phosphoenolpyruvate and carbon dioxide (bicarbonate) to form oxalacetate and phosphate ion. Oxalacetate is reduced to malate with simultaneous oxidation of an equimolar amount of reduced nicotinamide adenine dinucleotide (NADH) to NAD; the reaction is catalyzed by malate dehydrogenase (MDH). This results in a decrease in absorbance at 340 nm that is directly proportional to CO₂ concentration in the sample.

REAGENT COMPOSITION

- When reconstituted according to the direction, CO₂ reagent contains: PEP 1.8 mM, Magnesium Sulfate 10 mM, NADH 0.40 mM, MDH (porcine) 1.250 U/L, PEPC (microbial) 200 U/L, Sodium Oxamate 2.5 mM, Buffer, (pH 7.0). Non-reactive fillers and stabilizers with Sodium Azide 0.1% as preservative.
- CO₂ Standard contains 30 mmol/L of Sodium bicarbonate in an aqueous solution.

WARNINGS AND PRECAUTIONS

- CO₂ reagent is for "In Vitro Diagnostic Use."
- Normal precautions exercised in handling laboratory reagents should be followed.
- CO₂ reagent contains sodium azide that is toxic if ingested and which may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide accumulation.
- Use only CO₂-free water for reconstitution.
- Do not pipette by mouth to avoid CO₂ contamination from the expired air.

REAGENT PREPARATION

Reconstitute CO₂ Reagent with the volume of CO₂-free water, indicated on the vial label. Mix by inversion 5-6 times. Do not shake. CO₂ standard is ready to use.

Note: Avoid contamination of reagents with CO₂. Do not blow into pipette, since breath contains a high content of CO₂. Do not let bottles open unnecessarily, since CO₂ from air can contaminate reagent. Keep container tightly stoppered.

STORAGE AND STABILITY

- Store the unopened vial of CO₂ Reagent and CO₂ standard refrigerated (2-8°C). Reagent and standard are stable until the expiration date shown on the labels.
- Reconstituted reagents stored in the closed tightly stoppered bottle are stable for 1 day at room temperature (18-26°C) and 7 days refrigerated (2-8°C).

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Do not use the reagent if the absorbance of CO₂ reagent when measure against water in a 1-cm light path is less than 1 000 at 340 nm, or dry reagent exhibit caking due to possible moisture penetration or do not dissolve completely.

SPECIMEN COLLECTION AND STORAGE

1. It is recommended that specimen collection be carried out in accordance with NCCLS document M29-T.
2. Serum or heparinized plasma may be assayed. Ideally, venous blood should be collected and handled anaerobically.
3. Plasma and serum, after prompt separation from cells or clot, should be kept tightly stoppered.
4. CO₂ content of blood is stable for 1 hour when stored at 2-4°C under anaerobic conditions. 3.6 Alkalized serum can stand, without loss of CO₂, for as long as 4 hours in open cups.

INTERFERING SUBSTANCES

The major interference in this assay is CO₂ from air or from the breath of the analyst. Certain drugs and other substances are also known to influence blood CO₂ levels. See Young DS 8

MATERIALS PROVIDED

Carbon Dioxide Reagent and Carbon Dioxide Standard.

MATERIALS REQUIRED, BUT NOT PROVIDED

Spectrophotometer, capable of accurately measuring absorbances at 340 nm; Cuvettes with optical properties suitable for use at 340 nm; Pipetting devices for the accurate delivery of volumes required for this assay; Timer; Incubator or water bath, capable of maintaining 37°C.

PROCEDURE

1. Prepare CO₂ Reagent according to REAGENT PREPARATION.
2. Label tubes "Blank", "Standard", "Controls", "Patients", etc.
3. Pipette 1.0 mL Carbon Dioxide Reagent into each tube.
4. Incubate all tubes for 3 minutes at 37°C.
5. Set spectrophotometer wavelength at 340 nm, temperature to 37°C.
6. Pipette 5 µL (0.005 mL) of water, standard, and sample to the cuvette labeled "Blank", "Standard" and "Patients", respectively.
7. Mix gently by inversion and incubate for 5 minutes.
8. Read and record absorbance (Abs.) of all cuvettes at 340 nm.

PROCEDURE NOTES

If the samples is lipemic, icteric, or hemolyzed, a sample blank must be prepared by adding 5 µL of sample to 1.0 mL of saline. The absorbance of this blank is subtracted from the absorbance of the respective test and the corrected absorbance is then used in the calculations.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of normal and abnormal controls with known CO₂ values for each assay.

CALIBRATION

CO₂ Standard provided is standardized for using in manual procedure calculation. The concentration of CO₂ should be validated by comparison with commercially available standard. For automated instrument, refer to appropriate application specification.

CALCULATIONS

Determine CO₂ content of sample as follows:

$$\text{CO}_2 \text{ Content} = \frac{\text{Abs. blank} - \text{Abs. sample}}{\text{Abs. blank} - \text{Abs. standard}} \times \text{concentration of standard (mmol/L)}$$

Example:

Abs. sample = 0.95
Abs. standard = 0.82
Abs. blank = 1.66
Concentration of Standard = 30 mmol/L

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$$= 25,4 \text{ mmol/L}$$

PROCEDURAL LIMITATIONS

CO₂ reagent is linear to 40 mmol/L. If CO₂ content is greater than 40 mmol/L, dilute sample with an equal volume of CO₂-free isotonic saline and re-assay. Multiply the result by 2 to compensate for dilution. Lipemic, icteric, or hemolyzed sample will influence the result. Sample blank must be prepared. Please refer to PROCEDURE NOTES section for the preparation of sample blank.

EXPECTED VALUES

23-34 mmol/L

It is strongly recommended that each laboratory establish an expected range characteristic for the local population.

PERFORMANCE

1. **Linearity:** 0-40 mmol/L.

2. **Sensitivity:** Based on an instrument resolution of Abs = 0.001, this procedure has a sensitivity of 0.03 mmol/L.

3. **Comparison:**

Manual procedure: 40 serum samples with CO₂ range from 7.0-30.4 mmol/L were assayed. The study against similar method yielded the correlation coefficient of 0.987 and linear regression equation $y = 1.02x - 0.42$

Automated procedure: 70 serum samples with CO₂ range from 0.7-29.7 mmol/L were assayed. The study against similar method yielded the correlation coefficient of 0.99 and linear regression equation $y = 1.03x - 0.69$

4. Precision: The following studies were performed according to NCCLS guideline #EP5-T. Each control was subjected to 20 replicated assays for each study.

Within Run

	Control	Mean	S.D. (mmol/L)	C.V. (%)
Manual	I	34.0	0.6	1.8
Procedure	II	16	0.7	4.1
Automated	III	23.13	0.50	2.2
Procedure	IV	15.67	0.47	3.0

Run to Run

	Control	Mean	S.D. (mmol/L)	C.V. (%)
Manual	I	33.4	1.4	4.3
Procedure	II	15.7	0.7	4.4
Automated	III	23.29	2.00	8.6
Procedure	IV	14.79	1.60	

REFERENCES

1. Clinical Chemistry, LA Kaplan, AJ Pesce, Editors, CV Mosby Company, St. Louis (MO), p. 1056 (1984).
2. Hydrogen Ion Concentration in Body Fluids. In: Contarow and Trumper Clinical Biochemistry, 7th ed., AL Latner, Editor, Saunders, Philadelphia, p. 399 (1975).
3. Tietz NW, Pruden EL, Siggaard-Andersen O: Electrolytes, Blood Gases and Acid-Base Balance. In: Textbook of Clinical Chemistry, NW Tietz, Editor, Saunders, Philadelphia, p. 1188 (1986).
4. Natelson S: Routine use of ultramicro methods in the clinical laboratory. Am J Clin Pathol 21:1153, (1961).
5. Forrester RL, Wataji LJ, Silverman DA, Pierre KJ: Enzymatic method for the determination of CO₂ in serum. Clin Chem 22: 243, (1976).
6. Ibbott FA, LaGanga TS, Glin JB, Inkpen JA: Blood pH, CO₂ and

$$\text{CO}_2 \text{ (mmol/L)} = \frac{1,66 - 0,95}{1,66 - 0,82} \times 30 = 25,4 \text{ mmol/L}$$

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- O2. IN Clinical Chemistry- Principles and Technics, RJ Henry, DC Cannon, JW Winkelman, Editors, Harper & Row, Hagerstown (MD), p 755 (1974).
7. Gambino SR, Schreiber H: The measurement of CO2 content with the Auto Analyzer. A comparison with 3 standard methods and a description of a new method (alkalinization) for preventing loss of CO2 from open cups. Am J Clin Pathol 45:406, (1966).
8. Young DS: Effects of Drugs on clinical Laboratory Tests, 3rd ed., AACCPress, (1990).
9. Wolf PL: Interpretation of Biochemical Multitest Profiles. An Analysis of 100 Important Conditions. Masson Publishing, New York, p 33,56 (1977).

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METHOD: Immunoturbidimetric

FORMAT: Liquid

LINERITY: 2.3 % - 140 % (or value of highest calibrator)

SAMPLE: Plasma

KIT COMPOSITION:

R1 2 X 9.5 mL buffer

R2 1 X 6 mL antibody

Tests / Kit 120 (Hitachi 917)

(Can Vary Based On Analyzer Model)

INTENDED USE:

For the quantitative determination of Coagulation Factor XIII in human plasma.

REAGENT PREPARATION:

The reagents are ready to use and do not require reconstitution. Before use, gently invert Reagent 2 at least once a week.

STORAGE AND HANDLING:

All reagents should be stored refrigerated (2-8°C). Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for 18 months from the date of manufacture, as indicated by the expiration date on the package and bottle labels. Opened reagents can be used for one month if stored at 2-8°C.

PERFORMANCE:

Precision

The following within run precision results were obtained on a Roche/Hitachi analyzer with pooled human plasma:

	Sample 1	Sample 2	Sample 3
N	10	10	10
Mean	96.72 %	54.77 %	36.50 %
S.D.	0.6408	0.5945	0.2108
C.V.	0.66 %	1.09 %	0.58 %

Accuracy / Correlation

A comparison of the Factor XIII and another company's latex Factor XIII reagent was performed with the following results:

Linear Regression: $y = 0.9692x - 2.039$

$r = 0.9948$

$n = 50$

x = another company's latex Factor XIII

y = Factor XIII

Assay Range

2.3 % - 140 % (or value of highest calibrator)

INTERFERENCE:

Bilirubin C: No interference ≤ 22.0 mg/dL

Bilirubin F: No interference ≤ 19.7 mg/dL

Chyle (Formazine Turbidity): No interference $\leq 2,400$ FTU

Hemoglobin: No interference ≤ 450 mg/dL

Rheumatoid Factor: No interference ≤ 570 IU/mL

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Method: Immunoturbidimetric

Format: Liquid

Linearity: 2 - 80 µg/mL (or value of highest calibrator)

Sample: Plasma

KIT COMPOSITION

R1: 2 X 10.5 mL buffer

R2: 1 X 6 mL antibody

Tests / Kit: 120 (Hitachi 917)

(Can Vary Based On Analyzer Model)

Latex suspension, sera positive and negative controls,
glycine buffer 0.1 mol / L, pH 8.2;

INTENDED USE

For the quantitative determination of fibrin / fibrinogen
degradation products in human plasma.

PRINCIPLE

Latex particles with immobilized specific antibodies to
thrombolytic in the presence of degradation products of
fibrinogen / fibrin.

REAGENT PREPARATION

The reagents are ready to use and do not require
reconstitution. Before use, gently invert Reagent 2 at least
once a week.

STORAGE AND HANDLING

All reagents should be stored refrigerated (2-8°C). Return all
reagents to 2-8°C promptly after use. Unopened reagents
can be used for 18 months from the date of manufacture, as
indicated by the expiration date on the package and bottle
labels. Opened reagents can be used for one month if stored
at 2-8°C.

PERFORMANCE

Precision

The following within run precision results were obtained on a
Roche/Hitachi analyzer with pooled human plasma:

	Sample 1	Sample 2
N	20	20
Mean	5.185 µg/mL	18.515 µg/mL
S.D.	0.1694	0.2159
C.V.	3.27 %	1.17 %

Accuracy / Correlation

A comparison of the Plasma FDP reagent and another
company's Plasma FDP reagent was performed with the

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(Can Vary Based On Analyzer Model)

INTENDED USE:

The Ferritin assay is an *in vitro* diagnostic reagent for the quantitative determination of ferritin (an iron-storing protein) in human serum and plasma by immunoturbidimetric assay on the Roche / Hitachi 917 and other analyzers.

Measurements of ferritin aid in the diagnosis of diseases affecting iron overload and iron deficiency anemia.

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PERFORMANCE:

Precision

The precision was determined using packaged reagents, pooled human serum, and a Hitachi 717 analyzer

Precision - Within Run

Serum samples were assayed 21 times on the same day.

Sample I	Sample II	Sample III
N = 21	N = 21	N = 21
Mean = 14.9	Mean = 100.0	Mean = 431.1
Low = 13.8	Low = 98.7	Low = 427.6
High = 15.9	High = 101.2	High = 435.5
S.D. = 0.600	S.D. = 0.647	S.D. = 2.203
C.V. = 4.03%	C.V. = 0.65%	C.V. = 0.51%

Precision - Between Runs

Serum samples were assayed on 20 different days.

Sample I	Sample II	Sample III
N = 20	N = 20	N = 20
Mean = 16.5	Mean = 105.2	Mean = 428.7
Low = 15.4	Low = 102.7	Low = 421.0
High = 18.1	High = 107.6	High = 432.6
S.D. = 0.733	S.D. = 1.522	S.D. = 2.510
C.V. = 4.45%	C.V. = 1.45%	C.V. = 0.59%

Accuracy / Correlation

A comparison of the Ferritin and another company's Ferritin on the Hitachi 917 was performed with the following results. Statistics were calculated by linear regression. Samples used in the correlation study were serum samples from Asian patients.

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Fibrinogen Reagent Set Фибриноген

TYPE: Reagent

METHOD: Immunoturbidimetric
FORMAT: Liquid (Lyophilized Calibrator)

LINERITY:
100 - 900 mg/dL

or
81 - 729 mg/dL (using NIBSC international standard)

SAMPLE:
Plasma

KIT COMPOSITION:

Option 1

R1: 2 X 20 mL buffer
R2: 1 X 8 mL antibody
Calibrator: 1 x 1 mL (lyophilized)
Tests / Kit: 160 (Hitachi 717)
(Can Vary Based On Analyzer Model)

Option 2

Note: Large size kit, available by special order

R1: 4 X 50 mL buffer
R2: 1 X 40 mL antibody
(No calibrator)

Tests / Kit: 800 (Hitachi 717)
(Can Vary Based On Analyzer Model)

INTENDED USE:

For the quantitative determination of fibrinogen levels in disseminated intravascular coagulation (non-localized clotting within the blood vessels) and primary fibrinolysis (the dissolution of fibrinogen in a blood clot).
FOR IN VITRO DIAGNOSTIC USE.

REAGENT PREPARATION:

R1 and R2 are ready to use and do not require reconstitution. The calibrator is lyophilized and requires reconstitution and preparation of dilutions according to the package insert instructions.

STORAGE AND HANDLING

All reagents should be stored refrigerated (2-8°C). Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for one year from the date of manufacture, as indicated by the expiration date on the package and bottle labels. Opened R1 and R2 can be used for one month if stored at 2-8°C. Reconstituted fibrinogen calibrator and the diluted calibrator solutions can be used for 1 week if stored at 2-8°C.

PERFORMANCE:

Precision

Precision - Within Run

Sample I	Sample II	Sample III
N = 20	N = 20	N = 20
Mean = 142.4	Mean = 444.7	Mean = 739.6
Min = 136	Min = 429	Min = 721
Max = 152	Max = 456	Max = 763
C.V. = 2.96%	C.V. = 1.48%	C.V. = 1.43%

Precision - Between Runs

Sample I	Sample II	Sample III
N = 10	N = 10	N = 10
Mean = 167.4	Mean = 323.6	Mean = 554.1
Min = 165	Min = 316	Min = 544
Max = 170	Max = 329	Max = 571
C.V. = 0.94%	C.V. = 1.06%	C.V. = 1.50%

Accuracy / Correlation

A comparison of the Fibrinogen test and an Instar Fibrinogen Test Kit was

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INTENDED USE


For the quantitative determination of acid phosphatase in human serum

INTRODUCTION

Non-specific acid phosphatase activity is widely distributed throughout the living world. The enzyme secreted by human prostate gland has attracted most attention. Because of its clinical importance, extensive characterization and structural studies have now been carried out. Since acid phosphatase is also produced in other tissues, the prostatic isoenzyme must be distinguished from the non-prostatic for accurate diagnosis. Elevated levels of non-prostatic acid phosphatase have been observed in patients with Paget's disease, hyperparathyroidism with skeletal involvement, and in cancers that have invaded the bones.

Numerous phosphate compounds have been proposed as substrates for measuring acid phosphatase activity, such as phenylphosphate, p-nitrophenylphosphate, and thymolphthalein phosphate. α -Naphthylphosphate was proposed by Babson et al as a specific substrate for prostatic acid phosphatase. However, Amador et al demonstrated that this compound could be hydrolyzed by enzymes derived from other tissues. Hillman⁵ proposed a method in 1971 that included diazotized 2-amino-5-chlorotoluene (Fast Red TR) that formed a diazo dye that absorbed strongly at 405nm. L-tartrate was used as a specific inhibitor of prostatic acid phosphatase to establish differentially the amount of prostatic isoenzyme. The above kinetic method is specific, fast, simple and can easily be adapted to automated instrumentation.

PRINCIPLE

The α -Naphthol released from the substrate α -Naphthylphosphate by acid phosphatase is coupled with Fast Red TR to produce a colored complex that absorbs light at 405 nm. The reaction can be quantitated photometrically because the coupling reaction is instantaneous.

L-Tartrate inhibits prostatic acid phosphatase but does not interfere with the reaction mechanism. Therefore, if testing is performed in the presence or absence of L-Tartrate, the difference between the results of the two assays is the level of prostatic acid phosphatase in the serum.

REAGENT COMPOSITION

1. Acid phosphatase reagent (Concentrations refer to reconstituted reagent) α -Naphthylphosphate 3mM, Sodium Citrate 60mM, pH 5.3 $\pm$ 0.1.
2. L-Tartrate Reagent (Concentrations refer to reconstituted reagent): Sodium L-Tartrate 2M, Citric Acid 70mM, Sodium Citrate 10mM, pH 5.3 $\pm$ 0.1.
3. Acetate Buffer 5M, pH 5.0.

WARNINGS AND PRECAUTIONS

Reagents are for *in vitro* diagnostic use only.

REAGENT PREPARATION

1. Reconstitute acid phosphatase reagent with distilled water stated on the label. Swirl to dissolve.
2. Reconstitute L-Tartrate Reagent with 5.0 ml distilled water. Warm reagent to aid in dissolution, if necessary.
3. Acetate buffer is ready to use.

REAGENT STORAGE AND STABILITY

1. Unopened vials are stable until stated expiration date on vial label when stored refrigerated (2 - 8°C).
2. The reconstituted acid phosphatase reagent is stable for one day at room temperature (15 - 30°C) and for seven days when stored refrigerated at 2 - 8°C.
3. The reconstituted L-Tartrate Reagent is stable refrigerated at 2-8°C until expiration date indicated on the vial label. If crystallization of component occurs, warm at moderate temperature (40 - 50°C) until dissolved.
4. Acetate Buffer solution is stable refrigerated at 2 - 8°C until expiration date indicated on the vial label.

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1. The reconstituted acid phosphatase reagent, without serum added, has an absorbance greater than 0.4 when measured at 405nm against water.
2. The L-Tartrate Reagent is precipitated. Apply heat (40 - 50°C) to redissolve reagent.

SPECIMEN COLLECTION AND STORAGE

1. Use only clear, unhemolyzed serum.
2. Serum must be separated from clot within two hours after collection.
3. Acid phosphatase activity is extremely labile at room temperature. Stabilization of the enzyme can only be achieved by acidifying with the acetate buffer provided. Add 20 µl (0.02 ml) of buffer per 1.0ml of serum then mix. Treated serum samples will remain stable for seven days when kept refrigerated at 2 - 8 °C.
4. Do not use plasma. Some anticoagulants inhibit acid phosphatase activity and/or cause turbidity.

INTERFERING SUBSTANCES

1. High levels of bilirubin (icteric samples) reportedly inhibit acid phosphatase activity determined by this procedure.
2. A number of drugs and substances affect acid phosphatase activity. Young et al 10 has published a comprehensive list.

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1. Test tubes/tack
2. Accurate pipetting devices
3. Distilled/deionized water
4. Timer
5. Spectrophotometer capable of reading at 405 nm
6. Temperature must be closely controlled during assay. A temperature controlled (37°C) spectrophotometer cuvette should be used.

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MANUAL PROCEDURE

Note: Stabilize acid phosphatase immediately after separation of serum from the clot by adding 20µl (0.02ml) of Acetate Buffer per 1.0 ml of serum. Mix and store in refrigerator until assay is ready to be performed.

A. TOTAL ACID PHOSPHATASE

1. Reconstitute reagent according to instructions.
2. Label tubes, "Control", "Patient", etc.
3. Pipette 1.0ml of reagent into all tubes.
4. Zero spectrophotometer with water at 405 nm. Set cuvette temperature at 37°C.
5. Add 100 µl (0.10 ml) of sample to respective tube and allow incubating at 37°C for five minutes.
6. After incubation, read and record absorbance every minute for five minutes to determine ΔA/Minute.
7. Repeat procedure for each sample.
8. Values (U/L) are obtained by multiplying the ΔA/Minute by the factor. See "Calculations".

B. NON-PROSTATIC ACID PHOSPHATASE

1. Add 1.0 ml of reagent to appropriately labeled tube
2. Add 10 µl (0.01 ml) of L-Tartrate Reagent and mix.
3. Zero spectrophotometer with water at 405 nm. Set cuvette temperature to 37°C.
4. Add 100 µl (0.10 ml) of sample, mix and incubate at 37°C for five minutes.
5. After incubation, read and record absorbance every minute for five minutes to determine ΔA/Minute.
6. Values (U/L) are obtained by multiplying ΔA/Minute by the factor. See "Calculations".

C. PROSTATIC ACID PHOSPHATASE

The value is obtained by subtracting the result of the non-prostatic acid phosphatase assay (B) from the total acid phosphatase assay (A).

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1. The integrity of the reaction should be monitored by use of a normal and abnormal control serum with known acid phosphatase values.

2. Acid phosphatase in control sera is more labile than in fresh sera. Add 20 µl (0.02 ml) of acetate buffer per 1.0 ml of water used to reconstitute the control sera.

CALCULATIONS

One International Unit is defined as the amount of enzyme catalyzes the transformation of one micromole of substrate per minute under defined conditions.

A. Total Acid Phosphatase Calculation.

$$\frac{\Delta A / \text{min} \times 10^6 \times 1,1}{12,9 \times 10^3 \times 1,0 \times 0,1} = \text{U/L} = \Delta A / \text{min} \times 853$$

B. Non-prostatic Acid Phosphatase Calculation:

$$\frac{\Delta A / \text{min} \times 10^6 \times 1,1}{12,9 \times 10^3 \times 1,0 \times 0,1} = \text{U/L} = \Delta A / \text{min} \times 860$$

Where:

10⁶ = Conversion of moles to millimoles

1.1 = Total reaction volume (total A.P.)

1.11 = Total reaction volume (non-Prost A.P.)

12.9 x 10³ = Molar absorptivity of o-Naphthol-Fast Complex at 405 nm.

1.0 = Light path in cm

0.1 = Sample volume (ml)

SAMPLE CALCULATIONS

AA/min. total acid phosphatase = 0.01

AA/min. Non-Prostatic acid phosphatase = 0.009

Total acid phosphatase: 0.01 x 853 = 8.5 U/L

Non-Prostatic acid phosphatase: 0.009 x 860 = 7.7 U/L

Prostatic Acid phosphatase: 8.5 - 7.7 = 0.8 U/L

LIMITATIONS

Samples with values above 60 U/L at 37°C should be diluted 1:9 with normal saline, re-run, and the final result multiplied by 10.

EXPECTED VALUES

Total Acid Phosphatase: 0 - 9 U/L

Prostatic Acid Phosphatase: 0 - 3 U/L

Values were taken from literature¹¹. These values are referred to adults, both Males and females. It is strongly recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

PERFORMANCE CHARACTERISTICS

1. **Linearity.** 60 U/L at 37°C.

2. **Comparison.** A study performed using the method with a commercial reagent with a similar formulation yielded the following. N = 22.

	Total	Prostatic
Correlation Coefficient	0.97	0.98
Regression Equation	y = 0.96x + 0.38	y = 0.97x - 0.23

3. Precision.

	Within Run (N=20)	Total Acid Phosphatase
Mean (U/L)	8.9	19
S.D.	0.7	1.4
C.V. %	7.8	7.5
	Run to Run (N = 15)	Total Acid Phosphatase
Mean (U/L)	10.2	20.2
S.D.	1.2	1.5
C.D. %	11.7	7.4

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INTENDED USE

For the quantitative determination of alkaline phosphatase in human serum.

INTRODUCTION

Distributed in almost every tissue of the body, serum alkaline phosphatase (ALP) levels are of interest in the diagnosis of hepatobiliary disorder and bone disease. 1 Most of the ALP in the normal adult serum is from the liver or biliary tract. 2 Normal alkaline phosphatase levels are age-dependent, and are elevated during periods of active bone growth. Moderate elevations of ALP (not involving the liver or bone) may be attributed to Hodgkin's disease, congestive heart failure, and abdominal bacterial infections. 3 Elevations also occur in the third trimester of pregnancy.

Alkaline phosphatase is determined by measuring the rate of hydrolysis of various phosphate esters. p-Nitrophenyl Phosphate is one such ester that was used as a substrate by Fujita in 1939. 4 Bowers and McComb further modified the procedure to a kinetic assay. 5 In 1974, the Committee on enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology adopted a modification of the above procedure as the recommended procedure. 6 The present method is a modification of the reference methods of the above committee and the American Association for Clinical Chemistry. 7

PRINCIPLE

The enzymatic sequence employed in the assay of Alkaline Phosphatase is as follows.

Alk. Phos.

p-Npp is colorless but p-Nitrophenol has strong absorbance at 405 nm. The rate of increased absorbance at 405 nm is proportional to the enzyme activity.

REAGENT COMPOSITION

After combining R1 and R2 as directed, the reagent contains:
p-Nitrophenyl Phosphate 17 mM
Magnesium Ions 4 mM
Buffer (pH 10.2 ± 0.2)
Activator and binder

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.

CAUTION: *In vitro* diagnostic reagents may be hazardous.

Handle in accordance with good laboratory procedures which dictate avoiding ingestion, and eye or skin contact.

2. Specimens should be considered infectious and handled appropriately.

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STORAGE AND STABILITY

Store reagent set at 2-8°C (refrigerated). The reagents are stable until the expiration date if stored as directed. The working reagent is stable for 14 days at 2-8°C (refrigerated). Protect from direct light. Avoid microbial contamination.

REAGENT DETERIORATION

The reagent should be discarded if:

1. Turbidity has occurred, turbidity may be a sign of contamination.
2. The working reagent has an absorbance against water greater than 0.8 at 405 nm.

SPECIMEN COLLECTION

Unhemolyzed serum is the preferred sample. Heparinized plasma may also be used. Oxalate, fluoride and EDTA inhibit alkaline phosphatase, so are unsuitable as anticoagulants. 8 Samples should be kept cold and assayed as soon as possible after collection. A timed routine for sample collection and analysis should be established in each laboratory because ALP levels in serum or plasma, or in reconstituted control serum, rise significantly when stored at 2° - 8°C or at room temperature.

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INTERFERING SUBSTANCES

EDTA, citrate, fluoride, and oxalate inhibit alkaline phosphatase. Young et al. gives a list of drugs and other substances, which may interfere with the determination of ALP activity 9

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1. Pipetting devices
2. Test tubes/rack
3. Timer
4. Spectrophotometer with a temperature controlled cuvette
5. Heating bath/block

GENERAL INSTRUCTIONS

The reagent for alkaline phosphatase is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

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MANUAL PROCEDURE

1. Prepare working reagent according to instructions.
 2. Pipette 1.0 ml of reagent into appropriate tubes and allow to equilibrate to 37°C.
 3. Zero spectrophotometer with water at 405 nm.
 4. Transfer 0.025 ml (25 µl) of sample to reagent. Mix well.
 5. After one (1) minute, measure the absorbance. Return tube to 37°C. Repeat readings every minute for the next two (2) min.*
 6. Calculate the average absorbance difference per minute ($\Delta \text{Abs./min.}$)
 7. The $\Delta \text{Abs./min.}$ multiplied by the factor 2187 (see Calculations) will yield results in IU/L.
 8. Samples with values above 800 IU/L should be diluted 1:1 with saline, re-assayed and the results multiplied by two (2).
- *Note: If the spectrophotometer being used is equipped with a temperature controlled cuvette. The reaction mixture may be left in the cuvette while the absorbance readings are taken.

ALTERNATE VOLUMES

If the spectrophotometer being used can accurately read the final volume of 0.50 ml or smaller, follow the "ALTERNATE PROCEDURE".
Unit: One unit is the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under the specified conditions.

$$\text{IU/L} = \frac{\Delta \text{Abs. / min.} \times 1000 \times \text{TV}}{\square \times \text{LP} \times \text{SV}} = \frac{\Delta \text{Abs./min.} \times 1000 \times 1.025}{18.75 \times 1 \times 0.025} = \Delta \text{Abs./min.} \times 2187$$

Where,

$\Delta \text{Abs./min.}$ = Absorbance change
1000 = Conversion of IU/ml to IU/L
TV = Total reaction volume (1.025 ml)
 $\square$ = Millimolar absorptivity of p-Nitrophenol 18.75
LP = Light path in centimeters (1.0 cm)
SV = Sample volume (0.025 ml)
Example: If the $\Delta \text{Abs./min.}$ = 0.007 then $0.007 \times 2187 = 15.3 \text{ IU/L}$.
NOTE: If test parameters are altered the factor has to be recalculated using the above formula.
SI UNITS: To convert to SI Units (nkat/L) multiply IU/L by 16.67.

ALTERNATE PROCEDURE

1. Prepare working reagent according to instructions.
2. Pipette 0.50 ml (500 µL) of reagent into appropriate tubes and allow to equilibrate to 37°C.
3. Zero spectrophotometer with water at 405 nm.
4. Add 0.010 ml (10 µL) of sample to reagent, mix well.
5. After one (1) minute, measure the absorbance. Return tube to 37°C. Repeat readings every minute for the next two (2) min.*
6. Calculate the average absorbance difference per minute ($\Delta \text{Abs./min.}$)
7. The $\Delta \text{Abs./min.}$ multiplied by the factor 2720 (see Calculations) will yield results in IU/L.

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8. Samples with values above 800 IU/L should be diluted 1:1 with saline, re-assayed and the results multiplied by two (2).

*Note: If the spectrophotometer being used is equipped with a temperature controlled cuvette. The reaction mixture may be left in the cuvette while the absorbance readings are taken.

CALCULATIONS

Unit: One unit is the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under the specified conditions.

$$\Delta \text{Abs}/\text{min} \times 1000 \times 0.510 = \Delta \text{Abs}/\text{min} \times 2720$$

$$18.75 \times 1 \times 0.010$$

Where:

$\Delta \text{Abs}/\text{min}$ = Absorbance Change
1000 = Conversion of IU/ml to IU/L
TV = Total reaction volume (0.510 ml)
1 cm = Light path in centimeters
18.75 = Millimolar absorptivity of p-Nitrophenol
SV = Sample volume (0.010 ml)
Example: A1 = 0.58, A2 = 0.60
NOTE: (0.60 - 0.58) = 0.04 x 2720 = 108.8 IU/L
Then: If test parameters are altered the factor has to be recalculated using the above formula.
SI UNITS: To convert to SI Units (kat/L) multiply IU/L by 16.67.

LIMITATIONS

This methodology measures total alkaline phosphatase irrespective of tissue or organ of origin. Further tests may be necessary to assist in differential diagnosis.

QUALITY CONTROL

It is recommended that controls be included in each set of assays. Commercially available control material with established alkaline phosphatase values may be used for quality control. The assigned value of the control material must be confirmed, by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate reagent deterioration, instrument malfunction, or procedural errors.

EXPECTED VALUES

Adults 25 - 90 IU/L at 37°C. Children have a higher normal value. It is strongly suggested that each laboratory establish its own normal range 10

PERFORMANCE CHARACTERISTICS

1. Linearity: 900 IU/L.
2. Comparison: A group of 94 sera ranging in alkaline phosphatase values from 26.5 - 688 IU/L was assayed by this method and a similar commercially available reagent. Comparison of the results yielded a correlation coefficient of 0.997 and the regression equation was $y = 1.290x - 5.89$. (Comparison studies were performed according to NCCLS Tentative Guidelines, Ep9-T.)
3. Precision studies:

Within Run			Run-to-Run		
Mean IU/L	S.D.	C.V. (%)	Mean IU/L	S.D.	C.V. (%)
46.7	3.9	8.3	46.4	2.8	6.4
157.8	3.3	2.1	157.0	3.1	1.9

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Spectrophotometric Method for Measuring the Activity of Serum Alkaline Phosphatase. *Clin. Chem.* 12:70 (1966).

6. The Committee on Enzymes of The Scandinavian Society for Clinical Chemistry and Clinical Physiology: Recommended Methods for the Determination of Four enzymes in blood. *Scand. J. Clin. Lab. Invest.* 32:29 (1974).

7. Tietz, N.W., (ed): Study Group on Alkaline Phosphatase. A Reference Method for measurement of alkaline phosphatase activity in Human Serum. *Clin. Chem.* 29:751 (1983).

8. Demetrius, J.A., et al.: *Enzymes in Clinical Chemistry Principles and Techniques*, 2nd Ed. Hagerstown (MD), Harper and Row 927 (1974).

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Phosphorus Reagent Set Фосфор

INTENDED USE

Inorganic phosphorus (UV method) reagent is for the quantitative determination of inorganic phosphorus in human serum.

INTRODUCTION

The majority of the body's phosphorus is found in the bone as hydroxyapatite. The remaining phosphate is present as inorganic phosphate and phosphate esters. Phosphorus is involved in the intermediary metabolism of carbohydrates and is a component of other physiologically important substances. Thus, increased serum phosphorus may occur in hypervitaminosis, hypoparathyroidism, and renal failure. Reduced serum phosphorus levels are seen in rickets (vitamin D deficiency), hyperparathyroidism, and Fanconi's syndrome¹. The determination of inorganic phosphorus has been based on the reaction of molybdate with phosphate to produce the phosphomolybdenum blue complex, which is measured photometrically. However, many of the components in these reagents are unstable and had to be stored separately². Unreduced phosphomolybdate complex is measured directly in UV range 340 nm in the present method.

PRINCIPLE

Inorganic Phosphorus + H_2SO_4 + Ammonium Molybdate → Unreduced Phosphomolybdate Complex.
Inorganic phosphorus reacts with ammonium molybdate in an acid medium to form a phosphomolybdate complex, which absorbs light at 340 nm. The absorbance at this wavelength is directly proportional to the amount of inorganic phosphorus present in the sample.

REAGENTS

1. Inorganic Phosphorus Reagent: Ammonium Molybdate 0.4 mM, Sulfuric Acid 210 mM with surfactant.
2. Inorganic Phosphorus Standard: (5.0 mg/dl) Potassium Phosphate in dilute acid with an albumin base.

PRECAUTIONS

1. The reagents are for "in vitro" diagnostic use only.
2. Do not pipette by mouth. Avoid contact of reagents with skin, eyes and clothing.

REAGENT PREPARATION

Reagent comes in a ready to use form.

REAGENT STORAGE

Store reagent and standard at refrigerator temperature (2 - 8°C).

REAGENT DETERIORATION

1. Do not use if:
1. Reagent without serum added has an absorbance greater than in dilute acid with an albumin base.
2. The reagent fails to recover stated control values.

SPECIMEN COLLECTION AND STORAGE

1. Use only clear, unhemolyzed serum, separated from the erythrocytes as soon as possible. Erythrocytes contain organic phosphates that can hydrolyze on standing or can be enzymatically cleaved by phosphatases. Inorganic phosphates can then leak through the cell walls, increasing the concentration.
2. Once the serum has been separated, the phosphate content will not change for at least a week when stored in the refrigerator (2-8°C).

INTERFERENCES

For a comprehensive list of substances that interfere with the Measurements of Inorganic Phosphorus see Young, et al.

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Phosphorus Reagent Set Фосфор

1. Inorganic Phosphorus reagent
2. Inorganic Phosphorus standard

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1. Accurate pipetting devices.
2. Test tubes/rack.
3. Timing device.
4. Spectrophotometer with ability to read at 340 nm.

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PROCEDURE (MANUAL)

1. Label test tubes Blank, Standard, Control, Patient, etc.
2. Pipette 1.0 ml of reagent into each tube. Allow to come to room temperature (25°C).
3. Add 0.02 ml (20 µl) sample to respective tubes, mix, and allow to stand for five (5) minutes at room temperature.
4. Zero spectrophotometer with distilled water at 340 nm.
5. Read and record absorbencies of all tubes.

HTI - MULTI PURPOSE CALIBRATOR MAY BE USED TO REPLACE STANDARD.

PROCEDURE NOTES

1. For spectrophotometers requiring a larger total volume for accurate reading, a 3.0 ml reagent to 0.1 ml (100 µl) sample ratio may be used.
2. Lipemic and icteric samples require a serum blank. For maximum accuracy a serum blank should be run with each sample.
- a. Add sample to Saline solution.
- b. Zero spectrophotometer at 340 nm with Saline solution.
- c. Read and record absorbencies of serum blanks.
- d. Subtract absorbencies from test absorbencies.
3. Samples with values exceeding 12.0 mg/dl must be diluted 1:1 with saline, re-run, and result multiplied by two (2).

CALCULATIONS

Abs = Absorbance

$\text{Abs. of Unknown} - \text{Abs. of reagent Blank} \times \text{Cons. Of Stand} = \text{in org. Phosph}$

$\text{Abs. of Standard} - \text{Abs. of reagent Blank}$
= mg/dl

Example:

Abs. of reagent Blank = 0.536
Abs. of Unknown = 0.918
Abs. of Standard = 1.012
Cons. of Standard = 5 mg/dl

$$\frac{0.918 - 0.536}{1.012 - 0.536} \times 5 = \frac{0.382}{0.476} \times 5 = 4.0$$

To obtain results in SI units (mmol/L), multiply the result in mg/dl by the factor 0.323

Example: 4.0 mg/dl $\times 0.323 = 1.296$ mmol/L

LIMITATIONS

Most commonly employed detergent and disposable wipes used in the laboratory contain phosphates, and the use of improperly rinsed glassware may result in elevated inorganic phosphorus values.

CALIBRATION

It is not necessary to determine a standard curve since the reaction is linear in a range up to 12 mg/dl. However, a reagent blank and standard should be employed with each set of unknown assayed.

QUALITY CONTROL

It is recommended that controls be included in each set of assays.

Commercially available control material with established inorganic phosphorus values may be routinely used for quality control. The

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assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction or procedural errors.

EXPECTED VALUES

Adults: 2.5 - 4.8 mg/dl

This range should serve only as a guideline. It is recommended that each laboratory establish its own range of expected values since differences exist between instruments, laboratories and local populations.

PERFORMANCE

1. **Linearity:** 12 mg/dl

2. **Sensitivity:** Based on an instrument resolution of A - 0.001, the present procedure has a sensitivity of 0.01 mg/dl.

3. **Comparison:** A comparison study performed between this method and one based on the same methodology yielded correlation coefficient of 0.99 with a regression equation of $y = 1.01x - 0.06$.

4. **Precision:**

Day-to-Day Precision: Two commercial control sera were assayed for a period of thirty (30) days and the following day to day precision was obtained.

Day-to-Day (N = 21)

Mean (mg/dl)	S.D.	C.V. (%)
3.2	0.2	6.6
7.2	0.3	4.1

Within Run Precision: Two commercial control sera were assayed twenty (20) times and the following within run precision was obtained in org. Phos.

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Fructosamine Reagent Set Фруктозамин

INTENDED USE

This reagent is intended for the quantitative determination of fructosamine in human serum. For *in vitro* diagnostic use only.

CLINICAL SIGNIFICANCE

The determination of fructosamine is most commonly performed for the evaluation of glycemic control in diabetes. Fructosamine values provide an indication of glucose levels over the preceding 2-3 weeks. A higher fructosamine value indicates poorer glycemic control (1,2).

TEST SUMMARY

Glycated proteins are formed by a non-enzymatic reaction between glucose and protein in which unstable Schiff bases are formed, followed by an Amadori conversion to form stable ketamines. These glycated proteins include glycohemoglobin, glycoalbumin and glycated total protein.

Fructosamine is a term that has come into acceptance and refers to both glycoalbumin and glycated total protein. As the average life span of these proteins is about 2-3 weeks, the level of fructosamine provides a reflection of the average glucose concentration over that time (5).

Fructosamine and glycohemoglobin are both used to monitor diabetic control. However, each assay provides information for a specific time frame that is related to the analyte being measured.

Since the life span of hemoglobin is closer to 6-8 weeks, glycohemoglobin measurements reflect the average glucose concentration over this longer period of time (5).

Therefore, in comparison to glycohemoglobin determinations, fructosamine provides an index of intermediate-term diabetic control as opposed to the longer term for glycohemoglobin. Also, because of the shorter life span of the glycated albumin and total proteins, fructosamine measurements are more sensitive to changes in diabetic control. This provides a means to alert the physician to improvement, or deterioration in control much earlier than glycohemoglobin determinations (6).

There have been several methods developed for the determination of fructosamine. These methods include phenylhydrazine, furosine, affinity chromatography and several colorimetric procedures. A procedure using furosine and HPLC is accepted as the reference method; however, a colorimetric procedure using nitroblue tetrazolium (NBT) has gained popularity due to its speed, reproducibility and ease of automation (8). The reagent presented here is a modification of the commonly used NBT method. The fructosamine reagent set is based on the ability of ketamines to reduce NBT to a formazan dye under alkaline conditions.

The rate of formazan formation, measured at 550 nm, is directly proportional to the fructosamine concentration.

REAGENTS

Fructosamine Buffer: Carbonate buffer 100mM, pH 10.35 ± 0.1, sodium azide 0.1%.

Fructosamine Reagent: Nitroblue tetrazolium (NBT) 0.57mM, surfactant, non-reactive stabilizers and fillers.

Fructosamine Calibrator: Pooled human serum containing buffers, stabilizers and fillers.

REAGENT PREPARATION

Reconstitute the Fructosamine Reagent with the amount of Fructosamine Buffer specified on the vial label. Swirl gently to dissolve.

The Fructosamine Calibrator is supplied as a liquid stable, serum based product. It is ready to use upon opening.

REAGENT STORAGE AND STABILITY

1. Unreconstituted buffer, reagents and calibrators are stable until the expiration date on the kit label when stored at 2-8°C.

2. Upon reconstitution, the reagent should be stored at 2-8°C for best results.

3. Reconstituted reagent is stable for 7 days if stored at room temperature (15-25°C) or 30 days if stored refrigerated (2-8°C).

4. After opening, calibrator is stable for 30 days stored at 2-8°C.

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PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Do not use reagents past their expiration date stated on each reagent container label.
3. Do not pipette by mouth. Avoid ingestion and contact with skin.
4. Reagents in this kit contain sodium azide (0.1%) as a preservative. Sodium azide may form explosive compounds in metal drain lines. When disposing of reagents through plumbing fixtures, flush with copious amounts of water.
5. All specimens, calibrators and controls should be handled as potentially infectious, using safe laboratory procedures. (NCCLS M29-T2)9
6. Human serum was used in the manufacture of the calibrator. Each donor unit was tested and found negative, or non-reactive for HbsAg, HCV and HIV.

SPECIMEN COLLECTION AND STORAGE

1. Human serum, separated from the cells as quickly as possible, is the specimen of choice.
2. Collect specimens per NCCLS document H4-A3.10
3. Avoid hemolysis or contamination of the sample with hemoglobin as glycated hemoglobin will react in the same manner as fructosamine.
4. Serum specimens are stable for one week if stored at 2-8°C. Storage at -20°C is not recommended. 11

INTERFERENCES

1. All interference studies were performed according to the procedures recommended in NCCLS guideline No. EP7-P for interference testing in clinical chemistry. 12
2. Bilirubin to 20 mg/dl has been demonstrated to have a negligible effect (<5%) on fructosamine results using this method.
3. Hemoglobin to 200 mg/dl has been demonstrated to have a negligible effect (<5%) on fructosamine results using this method.
4. Glucose to 600 mg/dl has been demonstrated to have a negligible effect on fructosamine results using this method.
5. See Young, et al for other interfering substances. 13

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1. Pipetting devices.
2. Test tubes/rack.
3. Timing device.
4. Heating Block.
5. Spectrophotometer capable of reading at 550nm.
6. Fructosamine Controls, catalog number F7546-CTL.

TEST PROCEDURE (AUTOMATED-GENERAL)

Wavelength: 550nm
Assay Type: Kinetic with Standard
Sample/Reagent Ratio: 1:21
Reaction Direction: Increasing
Temperature: 37°C
Lag Time: 600 seconds
Read Time: 300 seconds
Reference Range: 1.61-2.68 mmol/L

TEST PROCEDURE (MANUAL)

1. Label tubes "Standard", "Control", "Sample", etc.
2. Pipette 1.0 ml reagent into all tubes and prewarm at 37°C for five minutes.
3. Add 0.05ml (50ul) of specimen to each respective tube at timed intervals.

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Fructosamine Reagent Set Фруктозамин

4. After exactly ten (10) minutes at 37°C read the absorbance of each tube at 550nm (A1). Return each tube to 37°C.
5. After exactly five (5) more minutes at 37°C, read tubes again at 550nm (A2).
6. To determine results see "Calculations".

LIMITATIONS

1. The procedure described is linear to 8.5 mmol/L. Samples with values exceeding 8.5 mmol/L should be diluted 1:1 with saline, re-assayed, and the result multiplied by two.
2. Hemoglobin greater than 200 mg/dl may give falsely elevated results.

CALCULATIONS

A = Absorbance

A2 Sample - A1 Sample x Conc. of Calibrator = Fructosamine

A2 Calibrator - A1 Calibrator in Sample

Example:

If A1 Sample = 0.100 and A2 Sample = 0.600,
A1 Calibrator = 0.100 and A2 Calibrator = 0.400,
and Concentration of Calibrator = 3.0 mmol/L, then:

$$\frac{0.600 - 0.100}{0.400 - 0.100} \times 3.0 \text{ mmol/L} = 5.0 \text{ mmol/L}$$

CALIBRATION

This assay requires the use of a fructosamine calibrator. The calibrator included in the kit is recommended.

QUALITY CONTROL

Serum controls with known normal and abnormal fructosamine values should be run routinely to monitor the validity of the reaction. These controls should be run at least with every working shift in which fructosamine determinations are performed. It is strongly recommended that each laboratory establish their own frequency of control determination.

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REAGENT PERFORMANCE

Assay Range: 1.0 - 8.5 mmol/L. Assay ranges for other analyzers may vary and must be validated.

Correlation: Results obtained with this reagent (y), in 45 samples ranging in fructosamine from 1.17-5.94 mmol/L, were compared with those obtained in the same samples using a reagent (x) based on the same methodology. The correlation coefficient was 0.988 and the regression equation was $y = 0.88x + 0.26$ (Std Err of Y Est = 0.19).

Precision: Precision studies were performed following a modification of the procedure contained in NCCLS document EP5-T2 15

Within Day (n=20)			Day to Day (n=20)		
Mean	S.D.	C.V. %	Mean	S.D.	C.V. %
1.97	0.04	2.0	1.91	0.06	3.1
5.57	0.10	1.8	5.72	0.14	2.4

Sensitivity: An investigation of the absorbance change per minute for ten replicates of two samples, with known concentrations of fructosamine, indicated that an absorbance change per minute of 0.042 was approximately equivalent to 1 mmol/L Fructosamine.

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6. Burtis, C.A., Ashwood, E.R., Tietz Fundamentals of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 370 (1996).
7. Burtis, C.A., Ashwood, E.R., Tietz Textbook of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 796 (1999).
8. Schleicher, E.D., Vogt, B.W., Clin. Chem. 36:136 (1990).
9. NCCLS document, "Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue", 2nd Ed. (1991).
10. NCCLS document, "Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture", 3rd Ed. (1991).
11. Burtis, C.A., Ashwood, E.R., Tietz Textbook of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 797 (1999).
12. NCCLS, "National Evaluation Protocols for Interference Testing", Evaluation Protocol Number 7, Vol. 4, No. 8, June 1984.
13. Young, D.S., et al, Clin. Chem. 21:1D (1975).
14. Burtis, C.A., Ashwood, E.R., Tietz Fundamentals of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 791 (1996).
15. NCCLS document, "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Ed. (1992).

- Burtis, C.A., Ashwood, E.R., Tietz Fundamentals of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 370 (1996).
- Burtis, C.A., Ashwood, E.R., Tietz Textbook of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 796 (1999).
- Schleicher, E.D., Vogt, B.W., Clin. Chem. 36:136 (1990).
- NCCLS document, "Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue", 2nd Ed. (1991).
- NCCLS document, "Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture", 3rd Ed. (1991).
- Burtis, C.A., Ashwood, E.R., Tietz Textbook of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 797 (1999).
- NCCLS, "National Evaluation Protocols for Interference Testing", Evaluation Protocol Number 7, Vol. 4, No. 8, June 1984.
- Young, D.S., et al, Clin. Chem. 21:1D (1975).
- Burtis, C.A., Ashwood, E.R., Tietz Fundamentals of Clinical Chemistry, Philadelphia (PA), W.B. Saunders, p. 791 (1996).
- NCCLS document, "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Ed. (1992).

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htimed.com

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Intended Use

For the quantitative determination of Chloride in serum.

Method History

In 1878 chloride was determined by precipitation of the anion as silver chloride and the measurement of the excess silver nitrate with standardized thiocyanate solution. Since then, many methods were introduced including mercurimetric titration using diphenylcarbazone as indicator. An electrometric method was introduced in the 1950's which was based on the exposure of a chloride solution to electrodes that emit silver ions in an acid medium. The time required for chloride precipitation as silver chloride is recorded and is directly proportional to chloride concentration. The present method was developed in the late 1950's and later adapted to the Technicon Auto Analyzer 6. This is a direct method and can be adapted to a variety of automated instruments.

Principle

The chloride ion displaces thiocyanate from non-ionized mercuric thiocyanate to form mercuric chloride and thiocyanate ions. The released thiocyanate ions react with ferric ions to form a color complex that absorbs light at 480nm. The intensity of the color produced is directly proportional to the chloride concentration.

Reagents

Chloride Reagent (Active Ingredients):
Mercuric Nitrate 0.105mM
Mercuric Thiocyanate 1.01mM
Ferric Nitrate 37.63mM
In dilute acid and methanol

Precautions

1. This reagent is for in vitro diagnostic use only.
2. This reagent is a TOXIC POISON as well as CAUSTIC. Avoid all contact with skin or clothing. Flush with water if contact occurs.
3. DO NOT PIPETTE BY MOUTH. Call physician if taken internally.

Reagent Preparation

Reagent comes in a ready to use form.

Reagent Storage

Store at room temperature protected from light.

Reagent Deterioration

Do not use if:
The reagent is a red-brown color and/or cloudy. The reagent should be a clear, pale-yellow solution.

Specimen Collection and Storage

1. Use serum that has been separated from the blood clot soon after drawing.
2. Grossly hemolyzed serum should not be used as it may create falsely decreased values.
3. Avoid contamination of blood with tissue fluid.
4. Store serum in tightly stoppered tubes.
5. Chloride is stable in serum for one day at room temperature and for three months frozen when stored tightly capped.

Interferences

1. Bromide and Fluoride can cause falsely elevated chloride values.
2. Other substances can influence chloride determination. For a comprehensive list see Young, D.S. et al.

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3. Lipemic and/or icteric serums do not interfere in the reaction

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1. Accurate pipetting devices.
2. Timer.
3. Test tubes/rack
4. Spectrophotometer able to read at 460-550nm.

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
Procedure (Manual)

1. Label test tubes "Blank", "Calibrator", "Sample", etc.
2. Add 1.0ml chloride reagent to each tube.
3. Add 0.01ml (10ul) of calibrator or sample to respective tubes, mix.
4. Incubate at room temperature for at least five minutes.
5. Set spectrophotometer at 480nm and zero with reagent blank. Wavelengths of 460-550nm may be used.
6. Read and record the absorbance reading of all vials.
7. See "CALCULATIONS" section to determine values. NOTE: Final color is stable for thirty minutes at room temperature.

Calibration

1. This reaction does not follow Beer's Law but has a linear range of 80-120 mEq/L on most instruments.
2. The linear range of the instrument to be used must be verified before actual testing. To verify a linear range, prepare a calibration curve by plotting the concentration of standards (ranging from 80-120mEq/L) versus their respective absorbance readings. Commercial standard/calibrator kits are recommended for this procedure.

3. Once the linear range has been established a single standard (included in the kit) can be used to standardize the reaction. Results are then calculated ratiometrically (see "CALCULATIONS").

Quality Control

Control serums with known concentrations of chloride should be used to monitor the integrity of the reaction.

Calculation

Abs = Absorbance
Abs of Unknown/ Abs. of Calibrator x Concentration of Calibrator
= Chloride (mEq/L)
Example: If a 100mEq/L calibrator has an absorbance reading of 0.550 and unknown has 0.450

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then:
0.450
— x 100 = 82 mEq/L
0.550

Limitations

1. Samples with chloride above 120 mEq/L should be diluted 1:1 with distilled water, re-run, and resulting value multiplied by two.
2. Do not touch pipette tips with fingers as chloride in skin perspiration may influence results.
3. Hydrochloric acid fumes may cause falsely elevated results.

Expected Values (R)

Serum 98-106 mEq/L

The above values are taken from literature and are intended as a guideline only. Each laboratory should establish its own range of expected values. Chloride levels in clinically healthy individuals vary by diet, sex, age, diurnal variation, physical activity, menstrual cycle, pregnancy and environmental factors.

Performance

1. Linearity: 80-120 mEq/L
2. Comparison: Studies done between this procedure and a similar procedure yielded a correlation coefficient of 0.991 with a regression equation of: $y = 1.01x + 0.11$.
3. Precision:

Within Run			Run to Run		
Mean	S.D.	C.V.%	Mean	S.D.	C.V.%
97	0.9	0.9	97	1.4	1.4
111	1.2	1.1	109	1.6	1.5

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INTENDED USE

For the quantitative determination of total cholesterol in serum.

INTRODUCTION

Cholesterol is a fatty substance found in blood, bile and brain tissue. It serves as a precursor to bile acids, steroids and vitamin D. The determination of serum cholesterol is a major aid in the diagnosis and classification of lipemias. Other conditions such as hepatic thyroid diseases influence cholesterol levels.²

Enzymatic methods have replaced older methodologies involving cholesterol esterase, oxidase, and trinders color system. Allain et al developed a total enzymatic technique in which hydrogen peroxide during the oxidation of cholesterol is used in conjunction with peroxidase, 4-aminoantipyrine and phenol to form a quinonimine dye.³ This reagent employs p-hydroxy benzene sulfonic acid (p-HBS), in place of phenol to produce a quinonimine dye with greater absorbance at 520 nm and a surfactant to facilitate the completion of reaction.

The results obtained from this replacement are equal to the results obtained based on the Abell-Kendall cholesterol reference method which is recommended by the CDC.

PRINCIPLE

The enzymatic reaction sequence employed in the assay of cholesterol is as follows:

Cholesterol Esters $\xrightarrow{\text{C. Esterase}}$ Cholesterol + Fatty Acids
Cholesterol + O₂ $\xrightarrow{\text{C. Oxidase}}$ Cholesterol-3-one + H₂O₂

2 H₂O₂ + 4-Aminoantipyrine
+ p-HBS

$\xrightarrow{\text{H. Peroxidase}}$ Quinonimine + 2 H₂O
(red dye)

Cholesterol esters are hydrolyzed to produce cholesterol. Hydrogen peroxide is then produced from the oxidation of cholesterol by cholesterol oxidase. In a coupled reaction catalyzed by peroxidase, quinonimine dye colored red is formed from 4-aminoantipyrine, p-HBS and hydrogen peroxide. The absorption at 520 nm of the solution of this dye is proportional to the concentration of cholesterol in the sample.

REAGENT COMPOSITION

Cholesterol (liquid) reagent set contains the following:

1. Cholesterol Reagent.
 - 4-Aminoantipyrine 0.6 mM, Sodium Cholate 8.0 mM, Cholesterol Esterase ≥ 150 U/L, Cholesterol Oxidase ≥ 150 U/L, Horseradish Peroxidase $\geq 1,200$ U/L, p-Hydroxy benzene sulfonate 20 mM, Buffer 125 mM, pH 6.8, non-reactive ingredients.
 2. Cholesterol Standard (Liquid): 200 mg/dl.
- This standard was made with materials traceable to Standard Reference Material available from the National Institute of Standard and Technology.

REAGENT PREPARATION

All reagents come in a ready-to-use form. No preparation is necessary.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
- CAUTION:** *In vitro* diagnostic reagents may be hazardous. Handle in accordance with good laboratory procedures which dictate avoiding ingestion, and eye or skin contact.
2. Specimens should be considered infectious and handled appropriately.

REAGENT STORAGE AND STABILITY

Store the reagent set at 2 - 8°C.

REAGENT DETERIORATION

1. The reagent should be discarded if:
1. Turbidity has occurred; turbidity may be a sign of contamination.
2. The reagent fails to meet linearity claims or fails to recover.

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control values in the stated range.

SPECIMEN COLLECTION

1. Test specimens should be serum and free from hemolysis
2. Cholesterol in serum is reported stable for seven (7) days at room temperature (18 - 25°C) and six (6) months when frozen and properly protected against evaporation.

INTERFERING SUBSTANCES

Anticoagulants such as fluoride and oxalate will result in false low values. The test is not influenced by hemoglobin values up to 200 mg/dl or by bilirubin levels up to 10 mg/dl. Interference from grossly icteric and heavily hemolyzed specimens is correctable by use of a serum blank.

MATERIALS PROVIDED

1. Cholesterol Liquid Reagent
2. Cholesterol Standard

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1. Spectrophotometer capable of measuring absorbances at 520 nm
2. Test tubes/rack
3. Accurate pipetting/measuring devices
4. Timer

GENERAL INSTRUCTIONS

The reagent for Cholesterol is intended for use either as an automated procedure on chemistry instruments or as a manual procedure on a suitable spectrophotometer.

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MANUAL PROCEDURE

1. Label test tubes: blank, standard, control, patient, etc.
2. Pipette 1.0 ml of reagent to all tubes and pre-warm at 37°C for at least 20 minutes.
3. Add 0.01 ml (10 µl) of sample to respective tubes, mix and return to 37°C.
4. Incubate all tubes at 37°C for ten (10) minutes.
5. Zero spectrophotometer with the reagent blank at 520 nm (Wavelength range: 500 - 550 nm).
6. Read and record absorbances of all tubes.

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NOTE

If the spectrophotometer being used requires a final volume greater than 1.0 ml for accurate reading, use 0.025 ml (25 µl) of sample to 3.0 ml of reagent. Perform the test as described above.

CALIBRATION

The procedures are calibrated with the standard solution, which is included with each series of tests. Its absorbance is used to calculate the results. It is recommended to establish a linearity curve up to 500 mg/dl with other available commercial standard solutions to verify the performance of instruments and reagents.

LIMITATIONS

1. The reagent is linear to 500 mg/dl.
1. Samples with values above 500 mg/dl should be diluted 1:1 with isotonic saline and re-run. Multiply final results by two (2).
2. Grossly lipemic serums require a "sample blank." Add 0.02 ml (20 µl) of sample to 2.5 ml saline, mix and read the absorbance against water. Subtract this value from the patient absorbance to obtain the corrected reading.

CALCULATIONS

(A = Absorbance)

$$\frac{A(\text{patient})}{A(\text{standard})} \times \text{Concentration of standard} = \text{Concentration of (mg/dl) patient (mg/dl)}$$

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A (patient) = 0.40, A (standard) = 0.32, Concentration of standard = 200 mg/dl

$0.40 \times 200 = 250 \text{ mg/dl}$
0.32

QUALITY CONTROL

1. It is recommended that high and low values of cholesterol controls be included in each set of assays. Commercially available control material with established cholesterol values may be used for quality control. The assigned value of the control material must be confirmed by the chosen application. Failure to obtain the proper range of values in the assay of control material may indicate either reagent deterioration, instrument malfunction, or procedural errors.
2. It is recommended that CDC or NIST certified reference material is included in the assay to assure the performance of the test.

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RISK CLASSIFICATION	TOTAL CHOLESTEROL IN BLOOD (MG/DL)
IN BLOOD (mg/dl)	
Desirable	< 200
Borderline high	200-239
High	≥ 240

PERFORMANCE CHARACTERISTICS

- Linearity:** 500 mg/dl.
- Sensitivity:** An absorbance change of 0.001 at 520 nm corresponds to 1 mg/dl under the stated condition of this assay system.
- Comparison:** A comparison between this procedure and one which is certified by the Center for Disease Control (CDC) and the National Cholesterol Education Program (NCEP) based on human samples assayed both by the Abell-Kendall's method and by the method being certified produced a regression equation of $y = 0.94x - 2.60$ ($N=47$) with a coefficient of correlation of 0.99.
- Precision:**

Within Run	S.D.	C.V. (%)
Mean (mg/dl)		
150	9.7	6.5
124	8.8	7.0

Run-to-Run	S.D.	C.V. (%)
Mean (mg/dl)		
142	7.4	5.2
114	6.1	5.3

5. **Specificity:** Cholesterol Oxidase is not totally specific for cholesterol. Other analogs of cholesterol (dihydrocholesterol, 7-dehydrocholesterol, 20 hydroxycholesterol, etc.) are also oxidized. These analogs do not normally occur in any appreciable amounts in serum.

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INTENDED USE


Cholinesterase reagent is used for *in vitro* diagnostic use in the quantitative kinetic determination of cholinesterase in human serum, plasma, or whole blood.

INTRODUCTION

Measurement of serum cholinesterase has been used to assess liver function and monitor excessive exposure to organophosphorus insecticides. It is also useful in predicting susceptibility to prolonged apnea after the administration of succinylcholine and to investigate the inheritance of variants of cholinesterase enzyme.^{1,2} There are many methods for the measurement of cholinesterase activity including manometric, titrimetric and photometric procedures. The colorimetric procedure based on the Ellman reaction is sensitive and simple. Therefore, this technique forms the basis of our reagent.

PRINCIPLE

Reactions involved in the cholinesterase assay are as follows:

Cholinesterase hydrolyzes Propionylthiocholine (PTC) to form Thiocholine which reacts with 5,5'-dithiobis-2-nitrobenzoic Acid (DTNB) to yield yellow 5-thio-2-nitrobenzoate with an absorbance maximum at 405 nm. Therefore, the rate of change in absorbance at 405 nm is directly proportional to cholinesterase activity.

REAGENT COMPOSITION

The Cholinesterase (PTC) Reagent, when reconstituted according to directions, has the following active ingredients. The approximate concentration of each component is as follows:

Propionylthiocholine 4mM, DTNB 0.4 mM, Buffer and non-reactive stabilizers and fillers

WARNING AND PRECAUTIONS

The Cholinesterase (PTC) Reagent is for *in vitro* diagnostic use.

Normal precautions exercised in handling laboratory reagents should be followed.

REAGENT STORAGE AND STABILITY

Store the dry reagent refrigerated (2 - 8°C). Reagent is stable until the expiration date shown on the label. Reconstituted reagent is stable for 6 hours at room temperature (15 - 30°C) or for 3 days at 2 - 8°C.

REAGENT PREPARATION

TOTAL CHOLINESTERASE ACTIVITY:

Reconstitute Cholinesterase (PTC) Reagent with volume of deionized water indicated on via label. After addition of water, stopper vial and immediately mix several times by inversion.

DIBUCAINE INHIBITOR:

Reconstitute one vial of Cholinesterase (PTC) Reagent with Dibucaine Solution, instead of deionized water.

FLUORIDE INHIBITOR:

Reconstitute one vial of Cholinesterase (PTC) Reagent with fluoride solution, instead of deionized water.

REAGENT DETERIORATION

The reagent should be discarded if:

1. Moisture has penetrated the vial and caking has occurred.
2. The reconstituted reagent has an absorbance against water greater than 1.2 at 405 nm.

SPECIMEN COLLECTION AND STORAGE

Cholinesterase activity is stable in undiluted serum for as long as 2 weeks at 2 - 8°C and up to 3 months at -20°C. Serum should be removed from clot promptly. EDTA does not inhibit cholinesterase activity.

Specimen preparation for measurement of cholinesterase activity in whole blood is as follows. Draw blood in tubes containing EDTA as anticoagulant. Mix whole blood thoroughly. An aliquot of blood is removed for hematocrit determination. Prepare blood hemolyzate by mixing 0.1 ml blood with 1.9 ml distilled water. Mix until hemolysis is complete. Centrifuge the remaining blood to obtain clear plasma. Determine the cholinesterase activities in plasma and hemolyzate following the instructions in "Procedure" section. Samples for hemolyzate cholinesterase activity are stable for 8 hours when stored at 2 - 8°C.

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INTERFERING SUBSTANCES

Hemolyzed serum samples should not be used. Certain drugs and other substances are also known to affect cholinesterase activity.

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1. Spectrophotometer, with temperature controlled compartment, capable of accurately measuring absorbance at 405 nm.
2. Test tubes with optical properties suitable for use at 405 nm.
3. Pipetting devices for the accurate delivery of volumes required for the assay.
4. Timer.
5. Dibucaine and Sodium Fluoride Solutions.

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MANUAL PROCEDURE

1. Reconstitute reagent according to instructions.
2. Pipette 1.0 mL of reagent into appropriate tubes and allow to equilibrate to 30°C or 37°C.
3. Zero spectrophotometer with water at 405 nm.
4. Add 10 µL of sample (serum, plasma or hemolysate). Mix well.
5. After 15 seconds, measure the absorbance (A_1). Return tube to 30°C or 37°C for another 30 seconds and measure another absorbance (A_2).
6. Calculate the ΔA per 30 seconds by subtracting A_1 from A_2 . Multiply by 2 to obtain ΔA per minute.
7. Calculate cholinesterase activity (U/L) of sample by multiplying $\Delta A/\text{min}$, times 7426.

ALTERNATE VOLUMES

The cholinesterase activity can also be determined by using 10 µL sample and 3 mL reagent. Different calibration factors should be used if sample and reagent volumes are different from the volumes required in the procedures.

DIBUCAINE AND FLUORIDE INHIBITION ASSAYS

Following the same steps given in the Procedure, determine cholinesterase activity in the samples using Dibucaine and Fluoride containing reagents prepared according to the instructions given in.

REAGENT PREPARATION: To determine percent of inhibition of activity, refer to **Calculations**.

CALCULATIONS (KINETIC)

$$\text{Cholinesterase Activity (U/L)} = \frac{\Delta A/\text{min} \times TV \times 1000}{13.6 \times LP \times SV}$$

$\Delta A/\text{min}$	= Absorbance change per min.
TV	= Total Volume (1.01 mL)
SV	= Sample Volume (0.01 mL)
13.6	= Millimolar absorptivity of DTNB
LP	= Lightpath (1 cm)
1000	= Conversion of units per mL to Units per Liter

$$\text{Cholinesterase Activity (U/L)} = \frac{\Delta A/\text{min} \times TV \times 1000}{13.6 \times 1.0 \times 0.01}$$

NOTE: Preparation of hemolysate involves 20-fold dilution of the sample; therefore, the Cholinesterase Activity in hemolysate calculated by the above formula should be multiplied by 20 to compensate for the dilution.

Example:

$$A_1 = 0.316 \quad A_2 = 0.491$$

$$A_2 - A_1 = 0.491 - 0.316 = 0.175$$

$$\Delta A/\text{min} = \Delta A/30 \text{ seconds} \times 2 = 0.350$$

Cholinesterase Activity (U/L) of sample = $0.350 \times 7426 = 2599$

DIBUCAINE AND FLUORIDE INHIBITION

Calculate percent of inhibition of Cholinesterase Activity (U/L) as follows:

$$\text{Dibucaine Inhibition of Cholinesterase Activity (\%)} = 100 - \frac{(CHE - D)}{CHE} \times 100$$

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Fluoride Inhibition of Cholinesterase Activity (%) = $100 - \frac{(CHE - F) \times 100}{CHE}$

Where:

CHE = Cholinesterase Activity in sample determined by using reagent containing neither Dibucaine nor Sodium Fluoride
CHE-D = Cholinesterase Activity in sample determined by using reagent containing 0.3 mMol/L Dibucaine
CHE-F = Cholinesterase Activity in sample determined by using reagent containing 40 mMol/L Sodium Fluoride
Erythrocyte Cholinesterase Activity
Erythrocyte Cholinesterase Activity (ECHE) is calculated from the results obtained for plasma Erythrocyte Cholinesterase Activity (PCHx), Hemolysate Cholinesterase (HCHE) and the hematocrit (HCT) value of the sample.

$$(HCHE) = (ECHE \times Hct) + (PCHx \times (1 - Hct))$$

$$ECHE = \frac{HCHE - (PCHx \times (1 - Hct))}{Hct}$$

* Hematocrit value expressed as decimal equivalent.

LIMITATIONS

1. Extremely lipemic or icteric serum should have blank correction performed.
2. This procedure does not include Dibucaine or Fluoride for resistance studies.

CALIBRATION

The procedure is calibrated by means of the millimolar absorptivity of 5-thio-2-nitrobenzoic acid, which is 13.6 at 405 nm. The linearity range of this reagent varies depending upon the sample to reagent ratio. The upper limit of linearity according to this procedure is 8,000 U/L with sample to reagent ratio of 1:100.

QUALITY CONTROL

It is recommended that high and low values of cholinesterase are included in each set of assays.
Commercially available control material with established cholinesterase values may be used for quality control.

Temperature Conversion Factor for Human Serum

Assay Temp.	Desired Temp	
25	1.00	1.20
30	0.83	1.00
37	0.65	0.77

EXPECTED VALUES

The expected range of serum cholinesterase activity at 30°C was determined to be as follows:

Serum	3100 -
7700 U/L	
Plasma	1700 -
4100 U/L	
Whole Blood	3300 -
5500 U/L	
Erythrocytes	4400 -
8200 U/L	

1. Linearity: 8000 U/L at 30°C.

2. Sensitivity: An absorbance change of 0.001/min. at 405 nm corresponds to 1.4 U/L of cholinesterase activity under the stated conditions of this assay system.

3. Comparison: A study performed between the present procedure and one commercial product on whole blood resulted in a coefficient of correlation of 0.96 with a regression of $y = 1.00x - 103$ ($n = 20$) and on serum/plasma resulted in a coefficient of 0.99 with a regression of $y = 0.96x + 63$ ($n = 47$).

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Intended use

Zinc test is an *in vitro* diagnostic assay for the quantitative determination of zinc in serum, plasma (not EDTA), cerebrospinal fluid (CSF) and urine. Only for *In Vitro* Diagnostics!

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Zinc test is a photometric assay showing the precision and accuracy of the atomic absorption spectrophotometry.

Zinc test is a kit of reagents that includes the sodium salt of 2-(5-Bromo-2-pyridylazo)-5-[N-propyl-N-sulfo-propyl]amino]phenol, (5-Br-PAPS) as the chromogenic component, and is used in a sensitive, colorimetric assay for zinc.¹ Wako Zinc test is employed to quantitate accurately the amount of zinc in a specimen such as serum, plasma, cerebrospinal fluid and urine.

Principle of the method

A solution of trichloroacetic acid is added to the specimen to precipitate proteins, leaving zinc in solution. After proteins are removed, Color Reagents A and B are added to an aliquot of the supernatant fluid. Zinc binds to 5-Br-PAPS, forming a reddish-violet chelate. The absorbance of this reddish-violet chelate is measured at 560 nm, and is directly proportional to the amount of zinc in the specimen. Although the color reagents employed can also be used for the determination of iron, copper, cobalt, and nickel, reaction with these other metals is suppressed by masking agents that are contained in the assay kit. 5-Br-PAPS chelates zinc with a 2 : 1 molar stoichiometry.

Reagents

Contents and storage conditions

R1:	Color Reagent A	Store at 2 - 25°C*
R2:	Color Reagent B	Store at 2 - 25°C*
R3:	Deproteinizing Reagent	Store at 2 - 25°C
CAL:	Standard Solution	Store at 2 - 25°C

* Store Color Reagents A and B under protection from light!

Ingredients

R1: 2 Bottles	Color Reagent A	60 mL
	Carbonate-Bicarbonate Buffer (pH 9.8)	0.2 mol/L
	5-Br-PAPS	0.07 mmol/L
	Sodium Citrate	170 mmol/L
	Dimethylglyoxime	4 mmol/L
R2: 1 Bottle	Color Reagent B	30 mL
	Salicylaldehyde	29 mmol/L
R3: 1 Bottle	Deproteinizing Reagent	30 mL
	Trichloroacetic Acid	7%
CAL: 1 Bottle	Standard Solution	10 mL
	Zinc	200 µg/dL

REAGENT PREPARATION and storage conditions

Working Solution:

Add 15 mL Color Reagent B (R2) to one bottle (60 mL) Color Reagent A (R1) and mix thoroughly. Other option: Mix Color Reagent A and Color Reagent B in a ratio of 4 : 1 (e.g., 20 mL of Color Reagent A is mixed with 5 mL of Color Reagent B). After mixing the Working Solution can be used for up to two days at room temperature (below 35°C), or for one week at 2 - 10°C.

After the preparation, the blank test value slightly increases over time. This has no effect on the measurement values.

All equipment must be free from zinc contamination. All reagent vials must be closed to avoid zinc contamination from airborne particles.

It is possible to use Color Reagents A and B individually. In this case, Color Reagent A must be added first.

SPECIMEN COLLECTION AND HANDLING

Serum, plasma (not EDTA), CSF and urine can be used as a specimen.

1. Take care to avoid contamination with chelators, such as EDTA, which cause artificially low zinc concentration values to be measured. Heparin, citrate and oxalate used as anticoagulants don't affect the measurement.

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High Technology, Inc

Zinc Reagent Set Цинк

2. In human blood, 75 - 85% of zinc exists in erythrocytes, 12 - 22% in plasma, and 3% in leukocytes. If hemolysis occurs, the amount of zinc that is measured in serum will be artificially high. In such cases, another sample of blood should be collected. Zinc levels can be determined in the serum, but also in plasma, CSF, and urine using a similar procedure.

3. Butyl rubber, which is often used as a stopper of the Control Serum bottle, contains some zinc. Dissolve the contents with „parafilm“ (or an equivalent) instead of butyl rubber as a stopper!

4. When sample of blood is collected, take care to avoid using apparatus** composed of rubber which contains zinc.

** ex. blood collecting tube, syringe.

MANUAL PROCEDURE

All equipment (Pipetting equipment, Aqua dest, Trichloroacetic Acid) must be free from zinc contamination. Wash sample and reagent containers 1 : 10 with Triton X-100 and water, immersed in dilute HCl or HNO₃ for 1 to 2 hours and rinse with distilled or deionized water. Invert and allow to drip dry.

Wavelength: 560 nm (Hg 546 nm).

Cuvettes: 1 cm light path.

Incubation Temperature: +15°-35°C.

Measure the absorbance of sample and standard against the reagent blank (RB). For each batch use reagent blank and one or two standards.

Details of procedure

Pipette into vials	Macro			Semi - Micro		
	Sample	Std.	RB	Sample	Std.	RB
Sample or control	0.5 mL			0.2 mL		
Standard		0.5 mL			0.2 mL	
Distilled or deionized water			0.5 mL			0.2 mL
Deproteinizing Reagent	0.5 mL	0.5 mL	0.5 mL	0.2 mL	0.2 mL	0.2 mL
Mix well. Centrifuge for 10 minutes at 3000 r.p.m.						
Supernatant	0.5 mL	0.5 mL	0.5 mL	0.2 mL	0.2 mL	0.2 mL
Working Solution	2.5 mL	2.5 mL	2.5 mL	1.0 mL	1.0 mL	1.0 mL
Mix well. Incubate for 10 minutes. Measure absorbance of the sample (A _s) and standard (A _{std}) against the reagent blank. The final color is stable for 1 hour.						

Limitation of the procedure

When zinc value exceeds 500 µg/dL, dilute sample 1 + 1 with saline or distilled water, repeat assay and multiply the result by 2.

Results

Calculations from the calibration curve The zinc concentration that corresponds to the measured absorbance of the sample (A_s) is determined by interpolation on a calibration curve.

$$\text{The zinc concentration in the sample } (\mu\text{g/dL}) = \frac{A_s}{A_{std}} \times 200.$$

Expected values

Serum and Plasma

Men: 72.6 - 127 µg/dL (11.1 - 19.5 mol/L)

Women: 70.0 - 114 µg/dL (10.7 - 17.5 mol/L)

Low levels may be encountered in women during pregnancy or menstruation.

Children: 63.8 - 110 µg/dL (9.8 - 16.8 mol/L)

New Born: 49.5 - 99.7 µg/dL (7.6 - 15.3 mol/L)

Urine: 300 - 800 µg/24 h

Since expected values are affected by age, sex, diet, geographical location and other factors.

Each laboratory should establish its own expected values for this procedure.

Performance characteristics

(1) Accuracy

The accuracy of the Zinc method was demonstrated by recovery

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High Technology, Inc

Zinc Reagent Set ЦИНК

studies in which known quantities of zinc were added to serum.

Recovery: 96 - 101%.

(2) Sensitivity

The theoretical sensitivity of this method expressed as absorptivity using a Hitachi 200 Spectrophotometer at 560 nm, is 1110 L/gm · cm for fresh reagent.

(3) Precision

The precision of the Wako Zinc method was established by analyzing a series of serum samples in replicate. CV% < 2%.

(4) Specificity

The specificity of the Wako Zinc method was determined by analysis of serum samples following the additions of potential interfering substances. Refer to tables 1 and 2.

Correlation^{2,3}

Specimen Serum

Correlation coefficient $r = 0.9228$ ($n = 82$)

Regression equation $y = 0.975x + 4.35$

y Wako Zinc 5-Br-PAPS Method (g/dL)

x Atomic Absorption Spectrophotometry (g/dL)

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Table 1: Serum Zinc Additive Study

Name	Additive		Serum with Additive Zinc Concentration (g/dL)
	Final Concentration	Ordinary Amount	
Bilirubin	None	-	71.3
	5 (mg/dL)	-	72.5
	10	-	70.5
	15	-	71.9
	20	-	72.5
	25	-	-
Hemoglobin	None	-	71.3
	200 (mg/dL)	-	77.6
	400	-	85.8
	600	-	92.1
	800	-	97.2
	1000	-	104.1
Anticoagulants	None	-	71.3
	2.0%	0.5%	72.5
	Citrate	0.01%	71.3
	Heparin	0.01%	72.6
	Oxalate	0.1%	1.9
	EDTA	0.1%	73.8

Table 2: Serum Zinc Additive Study (Effect of Metal Ion)

Name	Additive		Serum with Additive Zinc Concentration (g/dL)
	Final Concentration	Normal amount in Serum	
None	-	-	72.9
Cu ²⁺	100 (μg/dL)	100 (μg/dL)	72.9
	300	-	74.1
	500	-	75.3
	1000	-	77.6
Fe ³⁺	100 (μg/dL)	120 (μg/dL)	73.5
	300	-	74.1
	500	-	76.4
	1000	-	80.0
Fe ²⁺	100 (μg/dL)	120 (μg/dL)	73.1
	300	-	76.0
	500	-	77.6
	1000	-	81.7
Cd ²⁺	500 (μg/dL)	trace amount	76.4

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Warnings and precautions

- For in vitro diagnostic use.
- The usage and application of this test is reserved for professional use only. Please refer to respective national and local regulations and legislation.
- Not to be used internally in humans and animals.
- Do not use the reagents described above for any purpose other than described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Do not mix or use the reagents from one test unit with those of another test unit which has a different lot number!
- The reagent is designed to be used manually or on commercial available automated analyzers. Refer to the operating manual for a description of instrument operation and specifications.
- Clinical diagnosis must be determined with clinical symptoms and other test results by a physician.
- Store the reagents under the specified conditions. Do not use reagents past the expiration date stated on each reagent container label.
- After opening the reagents, it is recommended to use them immediately. When the opened reagents are stored, cap the bottles and keep them under the specified conditions.
- Since all specimens are potentially infectious, they should be handled with appropriate precaution. Refer to respective good laboratory practice protocols for preventing transmission of infection and handle samples in accordance with any other local or national regulations relating to the safe handling of such materials.
- If the reagents come in contact with the mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- When discarding the reagents, dispose of them according to local or national regulations.

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KIT COMPOSITION 1:

R1: 1 X 10 mL buffer
R2: 1 X 10 mL antibody
Tests / Kit: 100 (Hitachi 917)
(Can Vary Based On Analyzer Model)

KIT COMPOSITION 2:

R1: 2 X 10 mL buffer
R2: 2 X 10 mL antibody
Tests / Kit: 200 (Hitachi 917)
(Can Vary Based On Analyzer Model)

INTENDED USE:

For the quantitative determination of human cystatin C in serum, EDTA plasma, or lithium heparin plasma by immunoturbidimetric assay. Cystatin C measurements are used as an aid in the diagnosis and treatment of renal diseases. FOR IN VITRO DIAGNOSTIC USE.

REAGENT PREPARATION:

The reagents are ready to use and do not require reconstitution.

STORAGE AND HANDLING:

All reagents should be stored refrigerated (2-8°C). Return all reagents to 2-8°C promptly after use. Unopened reagents can be used for one year from the date of manufacture, as indicated by the expiration date on the package and bottle labels. Opened reagents can be used for one month if stored at 2-8°C.

PERFORMANCE:

The precision:

For the Cystatin C assay was determined using packaged reagents, control material, and a Roche / Hitachi 917 analyzer according to the CLSI EP5-A2 guideline.

	Sample 1 80	Sample 2 80	Sample 3 80	Sample 4 80
Mean (mg/L)	0.511	0.968	1.999	4.389

Within Run S.D.	0.006	0.007	0.013	0.030
Within Run C.V. %	1.094	0.712	0.640	0.690

Between Run S.D.	0.005	0.024	0.019	0.079
Between Run C.V. %	1.066	2.496	0.960	1.811

Between Day S.D.	0.005	0.017	0.014	0.023
Between Day C.V. %	0.928	1.776	0.707	0.525

Total S.D. 0.007	0.024	0.027	0.088	
Total C.V. %	1.421	2.462	1.353	2.008

Accuracy / Correlation:

Testing was performed on a Roche / Hitachi 917 analyzer according to the CLSI EP5-A2 guideline. A comparison of the Cystatin C and another company's Cystatin C assay was performed with the following results.

Linear Regression: $y = 1.0093x - 0.0411$

$r = 0.9983$

$n = 50$

x = another company's cystatin C assay

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
Summary and Explanation of Test

Ketosis is a common feature in acutely ill patients. In subjects suffering from starvation, acute alcohol abuse, or diabetes mellitus ketosis can result in severe life threatening metabolic acidosis. The presence and degree of ketosis can be determined by measuring blood levels of β -hydroxybutyrate. Ordinarily, β -hydroxybutyrate is the ketoacid present in the greatest amount in serum. It accounts for approximately 75% of the ketone bodies that also contain acetoacetate and acetone. 2,3,4 During periods of ketosis, β -hydroxybutyrate increases even more than the other two ketoacids, acetoacetate and acetone, and has been shown to be a better index of ketocidosis including the detection of subclinical ketosis. 5,6,7,8 In diabetics, the measurement of β -hydroxybutyrate as well as blood glucose is needed for the assessment of the severity of diabetic coma and is essential for the exclusion of hyperosmolar non-ketotic diabetic coma. Moreover, the insulin requirements are often based on the extent of the existing hyperketonemia shown by the blood levels of β -hydroxybutyrate. 9 A specific enzymatic assay for β -hydroxybutyrate is therefore extremely important in the assessment of ketosis.

Principle

Enzymatic quantitation of β -hydroxybutyrate by β -hydroxybutyrate dehydrogenase has been reported. In the High Technology Inc. method, β -hydroxybutyrate (D-3-hydroxybutyrate) in the presence of NAD gets converted to acetoacetate and NADH at a pH 8.5 by β -hydroxybutyrate dehydrogenase (D-3-hydroxybutyrate dehydrogenase). At this pH, the reaction is favored to the right. 12 The NADH produced is converted to color using INT and diaphorase.

Reagents

The 3-hydroxybutyrate reagent kit contains the following:

- 1 x 50ml R1 containing 3-hydroxybutyrate dehydrogenase and diaphorase enzymes.
- 1 x 8ml R2 containing NAD, INT and oxalate.
- 1 x 3ml standard containing 1mM sodium D-3-hydroxybutyrate

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Reagent Preparation

R1 and R2 reagents are supplied as ready to use liquids. To prepare a single liquid working reagent mix 10 parts R1 with 1.5 parts R2 (ex. 10ml R1 with 1.5 ml R2).

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Interferences

Lactate dehydrogenase and lactate have been shown to interfere with the assay. The incorporation of oxalic acid in this reagent eliminates this interference as reported. 12 No significant changes in values were observed when the following analytes were added to serum containing 0.5mM β -hydroxybutyrate.

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1. 10 ml Micropipettor
2. 1 ml Pipette or Dispenser
3. Glass or plastic Test Tubes to hold 1.0 ml
4. β -hydroxybutyrate controls
5. Temperature controlled incubator
6. Spectrophotometer capable of reading at 505nm.

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Sample volume 25 ul

Incubation 5 min.

1. Prepare working reagent according to "Reagent Preparation" instructions.
2. Label tubes, "Standard", "Control", "Patient", etc.
3. Pipette 1.0 ml of working reagent into required number of labeled test tubes. Incubate at 37°C for 5 minutes.
4. Add 25 ul Standard, controls and patients to the appropriate tests and incubate at 37°C for 5 minutes.
5. Zero spectrophotometer with a reagent blank.
6. Read and record absorbances of all tubes at 505 nm.

NOTES:

- a) Precise measurement of temperature, wavelength, and time are required to obtain accurate results.
- b) The test can also be performed at 25°C for 10 minutes using the same procedure as above.

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Quality Control

The reliability of test results should be monitored whenever patient samples are assayed using a standard and quality control materials analyzed in the same manner employed for the unknowns. We suggest the use of commercially available β -hydroxybutyrate controls with an assayed range, if controls do not fall into the assayed range, patient values from that run should not be reported. The run should be repeated, making sure that all mixing and handling instructions are strictly followed.

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Expected Values

The quantitation of β -hydroxybutyrate is important in cases of ketoacidosis. In studies of healthy individuals who had fasted for 12 hours before blood collection, the range of β -hydroxybutyrate was found to be from 0.02 mM (0.2 mg/dl) to 0.27 mM (2.81 mg/dl). 4.5 Other ranges have also been reported.

Performance

- Linearity: The β -hydroxybutyrate assay is linear to 4.5 mM (46.8 mg/dl).
- Comparison: A study of 57 samples using High Technology Inc. and the Sigma method gave the following correlation: $(y = \text{High Technology Inc.}, x = \text{Sigma})$ $y = 1.05x + 0.25$ mM, $N = 57$, $r = 0.9922$, $S_{xy} = 0.1055$ mM.
- Precision: Studies were conducted using two serum pools containing 0.25 mM (2.6 mg/dl) and 1.0 mM (10.4 mg/dl) β -hydroxybutyrate. The following results are averages of eighteen determinations:

Within day:

Level	Mean (mM)	Std. Dev. (mM)	% C.V.
1	0.29	0.005	1.7
2	1.09	0.015	1.4

Day to day:

Level	Mean	Std. Dev.	% C.V.
1	0.26	0.014	5.2
2	1.05	0.018	1.7

- Sensitivity: Concentrations of β -hydroxybutyrate of 0.18, 0.28, and 0.38 mM (1.8, 2.9, and 3.9 mg/dl) can be clearly distinguished at the 99% confidence limit.

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Intended Use

For *in vitro* diagnostic use only.
The HTI, Inc. Ammonia / Alcohol Control set is to be used for monitoring the accuracy and precision of various ammonia and / or ethanol assay methods and to validate quantitation of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The HTI, Inc. control material is supplied as a two level control set, 2 x 5 ml, as a ready-to-use liquid requiring no reconstitution or dilution. It is prepared in a aqueous base fortified with ethanol and reagent grade chemicals. Preservatives have been added to inhibit growth.

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Storage Stability

	Storage	Stability
Unopened	2-8°C	Refer to the label on each vial and on the package for the expiration date
Open Vial Stability	2-8°C	20 days

General Instructions for Use
Use the quality control material according to the directions accompanying the instrument or the assay procedure used. Treat the quality control material in the same manner as patient samples.

1. Gently invert the vial and swirl to assure contents are homogenous.
2. Remove the screw cap and gently remove the rubber stopper from the vial.
3. Perform assay.
4. Record the results according to your quality assurance program.

Expected Results

Refer to the expected results supplied for mean and range value assignments. See table below. Verify vial lot number with those on this assay sheet.

Assay	Unit	Level 1		Level 2	
		Mean	range	mean	range
Ammonia *	umol/L	48	35-61	244	195-293
Alcohol **	mg/dl	53	45-61	262	220-304
	mmol/L (SI)	11.5	9.7-13.3	56.8	47.7-65.9

Data generated using Hitachi 917 and Cobas Mira.

Data generated using Hitachi 717 and Cobas Mira.

The assay values and expected ranges are target values derived from inter laboratory data. The expected range values include variations of instrument and laboratory handling. The assay values were obtained using in-date HTI reagents available at the time of testing. Updates to the listed values may be made based upon additional data that becomes available or, if necessitated by a modification to a test method. The mean values established for your laboratory should fall within the ranges shown in Expected Results; however laboratory means may vary during the life of the control. Each laboratory should establish its own mean and precision parameters.

Limitations

The results obtained using the quality control material depends on several factors: erroneous results can occur from improper storage, inadequate mixing, or sample handling errors associated with instrument or assay procedures. Do not use the quality control material if there is visible evidence of microbial growth in the vial.

For more information about procedural limitations, refer to your instrument manual or assay product insert.

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Manufactured by
HTI MEDICAL 109
PRODUCTION RD
WALPOLE, MA 02081 U.S.A.

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Тип: Control

KIT COMPOSITION 1:

Level 1 Human Serum 1 x 1 mL, lyophilized
Level 2 Human Serum 1 x 1 mL, lyophilized

KIT COMPOSITION 2:

Level 1 Human Serum 2 x 1 mL, lyophilized
Level 2 Human Serum 2 x 1 mL, lyophilized

INTENDED USE:

The Apolipoprotein Control is intended for use as a consistent test sample of known concentration for monitoring assay conditions for the listed constituents. Assay values are provided for the specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and instrument.

CONTROL PREPARATION:

1. Allow controls to come to room temperature. Remove cap carefully.
2. Add 1.0 mL of deionized water. Let rest for 15 minutes.
3. Invert gently until dissolved.

STORAGE AND HANDLING:

Store controls at 2-8°C. After reconstitution, the controls are stable for 2 weeks at 2-8°C or 2 months if aliquoted and stored at -20°C. Unopened controls can be used up to the expiration date shown on the vial label and package insert.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Apolipoprotein Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

ANALYTE (mg/dL)	MEAN (mg/dL)	LEVEL 1 RANGE (mg/dL)	LEVEL 2 RANGE (mg/dL)	MEAN (mg/dL)	RANGE
Apolipoprotein AI	116	93 - 139	172	146	116 - 198
Apolipoprotein AII	27	22 - 32	42	36	27 - 48
Apolipoprotein B	53	42 - 64	105	89	53 - 121
Apolipoprotein CII	2.7	2.2 - 3.2	5.1	4.3	2.7 - 5.9
Apolipoprotein CIII	8.5	6.8 - 10.2	14.6	12.4	8.5 - 16.8
Apolipoprotein E	3.1	2.5 - 3.7	5.3	4.5	3.1 - 6.1

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Тип: Control

KIT COMPOSITION 1:

Level 1 Human Serum 1 x 2 mL, lyophilized
Level 2 Human Serum 1 x 2 mL, lyophilized

KIT COMPOSITION 2:

Level 1 Human Serum 2 x 2 mL, lyophilized
Level 2 Human Serum 2 x 2 mL, lyophilized

INTENDED USE:

The ASO/RF/CRP Control is intended to be used as a consistent test sample of known concentrations for monitoring the Anti-Streptolysin O (ASO), Rheumatoid Factor (RF), and C-Reactive Protein (CRP) immunoturbidimetric assays. FOR IN VITRO DIAGNOSTIC USE

CONTROL PREPARATION:

- Carefully and slowly remove rubber stopper and let air enter vacuum-sealed bottle.
- Carefully add exactly 2 mL purified water (at room temperature) to bottle. DO NOT MIX.
- Replace rubber stopper.
- Let bottle sit at room temperature for 10 minutes.
- Mix by slowly swirling bottle numerous times. DO NOT SHAKE or induce foaming.
- Continue to mix until all contents have dissolved completely.
- Before each subsequent use, be sure to mix bottle by slowly swirling.

STORAGE AND HANDLING:

Store lyophilized and reconstituted controls at 2-8°C. Return all controls to 2-8°C promptly after use. Unopened controls can be used for one year from the date of manufacturing as indicated by the expiration date on the package and bottle labels.

Reconstituted controls can be used for 2 weeks if stored at 2-8°C. Reconstituted controls can also be frozen immediately after reconstitution at -20°C for 2 months and then thawed ONE TIME ONLY.

EXPECTED VALUES:

The following values are given as an example. The expected values for the ASO/RF/CRP Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

ANALYTE		LEVEL 1		LEVEL 2	
MEAN	RANGE	MEAN	RANGE	MEAN	RANGE
ASO	258 IU/mL	232 - 284 IU/mL	427 IU/mL	384 - 470 IU/mL	
RF	32 IU/mL	22 - 42 IU/mL	82 IU/mL	66 - 98 IU/mL	
CRP	1.1 mg/dL	0.8 - 1.4 mg/dL	3.8 mg/dL	3.3 - 4.3 mg/dL	
CRP(3) Cal D	13.62 mg/L	10.89 - 16.34 mg/L	41.20 mg/L	32.96 - 49.44 mg/L	
CRP(3) Cal E	11.12 mg/L	8.89 - 13.34 mg/L	37.41 mg/L	29.93 - 44.89 mg/L	
CRP(3) Cal F	11.5 mg/L	9.2 - 13.8 mg/L	37.1 mg/L	29.7 - 44.6 mg/L	

ASO: Standardized against the WHO NIBSC Anti-Streptolysin-O standard material.

RF: Standardized against the WHO International Reference Preparation of Rheumatoid Arthritis Serum.

CRP, CRP(3): Standardized against CAP/BCR/IFCC Reference Preparation for Proteins in Human Serum (RPHS), lot 91/0619.

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Intended Use

For in vitro diagnostic use only

The Control Set is to be used for monitoring the accuracy and precision of various B-Hydroxybutyrate dehydrogenase assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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	Storage	Stability
Unopened	2-8°C	Refer to the label for the expiration date
Open Vial Stability	2-8°C	20 days

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Disposing of Materials

Dispose of hazardous or biologically contaminated materials according to your institution's practices. Discard all materials in a safe and acceptable manner that is in compliance with all county, state, and local requirements.

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High Technology, Inc

B-Hydroxybutyrate dehydrogenase Control Set

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High Technology, Inc

Glycosylated hemoglobin (HbA1c) Control Set

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INTENDED USE:

This product is for the purpose of monitoring accuracy and precision in the quantitative determination of human Hemoglobin A1c (HbA1c) in blood by automated immunoassay FOR IN VITRO DIAGNOSTIC USE.

CONTROL PREPARATION:

1. Reconstitute vials with 0.5 mL deionized water. Gently mix for 10 minutes. Observe for undissolved material.
2. The reconstituted controls may be dispensed in 0.1 mL aliquots, sealed tightly and frozen at -20°C.

STORAGE AND HANDLING:

1. Store controls at 2-8°C. Stable until expiration date if sealed tightly. PROTECT FROM LIGHT AND HEAT.
2. Reconstituted controls retain their assigned values for at least three months if frozen. If not frozen, the reconstituted controls are stable at least one month if stored at 2-8°C and sealed tightly.
3. Do not freeze and thaw more than once.
4. Do not store in a self-defrosting freezer.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Hemoglobin A1c Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

CONTROL LEVEL	% NGSP Units (RECOMMENDED)	mmol / mol Hb		
IFCC Units	MEAN	RANGE	MEAN	RANGE
Level 1	5.0 %	4.0 - 6.0 %	31	25 - 37
Level 2	11.9 %	10.4 - 13.4 %	107	93 - 121

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Intended Use

Controls for use in the quantitative determination of human glycohemoglobin (HbA1c) in blood by cation exchange resin.

Summary and Principle

Controls should be included each time patients are assayed for glycohemoglobin to verify that the assay has worked correctly. The mean value of the controls were obtained by assaying representative samples of the entire lot.

Reagents

The lyophilized glycohemoglobin controls are hemolysates prepared from packed human erythrocytes. The controls provide two levels of glycohemoglobin in the normal and elevated range. Stabilizers are added to maintain hemoglobin in the reduced state providing complete control of the glycohemoglobin procedure.

Reagent Preparation

Reconstitute vials with 1.0 ml deionized water. Gently mix for 10 minutes. Observe for undissolved material. The reconstituted controls should be dispensed in 0.1 ml aliquots, sealed tightly and frozen at -20°C.

Limitations

Things to look for that might cause inaccurate results are improper pipetting, inadequate mixing and poorly calibrated instruments.

Expected Values

See values listed on vial label. The assayed limits are to be used as a guide in determining the accuracy of the assay procedure. The assay results for the controls should fall within the stated expected range. If they do not, the assay should be repeated, checking closely for the factors mentioned in "Limitations".

Reagent Storage

1. Store reagent at 2-8°C. Stable until expiration date if sealed tightly.
PROTECT FROM LIGHT AND HEAT.
Reconstituted controls retain their assigned values for at least three months if frozen. If not frozen, the reconstituted controls are stable at least one month if stored at (2-8°C) and sealed tightly. Do not freeze and thaw more than once. Do not store in a self-defrosting freezer.

P

Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Although this product has been tested and found non-reactive for Hepatitis B Surface Antigen (HbsAg) and HIV, no known test can offer assurance that products derived from human blood will not transmit disease. Therefore all human serum products and patient specimens should be handled in the same manner as an infectious agent.
3. Do not pipette by mouth. Avoid contact with skin and mucous membranes.

Procedure

The lyophilized glycohemoglobin controls should be assayed in the same manner as blood specimens including the hemolysis procedure. Follow the directions that accompany the instrument and reagent kit used in the assay.

Materials Provided

1. Normal level control
2. Elevated level control

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Manufactured by:
HTI Medical
109 production Rd
Walpole, MA 2081

Произведено High Technology Inc.
109 Production Rd
Walpole, MA 02081 USA
Fax: 1 508 660-22-24
e-mail: medtech@htimed.com

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Intended Use

For *in vitro* diagnostic use only

The Control Set is to be used for monitoring the accuracy and precision of various Glucose-6-phosphatodehydrogenase (G6PD) assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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General Instructions for Use

Use the quality control material according to the directions accompanying the instrument or the assay procedure used. Treat the quality control material in the same manner as patient samples. Gently invert the vial and swirl to assure contents are homogenous. Remove the cap from the vial. Perform assay. Record the results according to your quality assurance program.

Expected Results

Refer to the expected results supplied for mean and range value assignments. Verify vial lot number with those on the assay values sheet. The assay values and expected ranges are target values derived from inter laboratory data. The expected range values include variations of instrument and laboratory handling. The assay values were obtained using reagents available at the time of testing. Updates to the listed values may be made based upon additional data that becomes available or, if necessitated by a modification to a test method. The mean values established for your laboratory should fall within the ranges shown in Expected Results, however laboratory means may vary during the life of the control. Each laboratory should establish its own mean and precision parameters.

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INTENDED USE:

The Liquid ITA Control is intended to be used as a consistent test sample of known concentration for monitoring the performance of immunoturbidimetric assays (ITA) for the listed constituents.
FOR IN VITRO DIAGNOSTIC USE.

CONTROL PREPARATION:

Controls are liquid and ready to use.

STORAGE AND HANDLING:

Store controls at 2-8°C. Return all controls to 2-8°C promptly after use. Unopened controls can be used for one year from the date of manufacturing as indicated by the expiration date on the package and bottle labels. Opened controls can be used for 1 month if stored at 2-8°C.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Liquid ITA Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

ANALYTE	LEVEL 1		LEVEL 2	
	MEAN (mg/dL)	RANGE (mg/dL)	MEAN (mg/dL)	RANGE (mg/dL)
RANGE (mg/dL)				
Alpha-1 Acid Glycoprotein	69	55 - 83	130	117 - 143
Alpha-1 Anti-Trypsin	115	92 - 138	224	202 - 246
C-Reactive Protein	1.1	0.8 - 1.4	3.8	3.3 - 4.3
Complement C3	98	78 - 118	243	219 - 267
Complement C4	19	15 - 23	45	40 - 50
Haptoglobin	96	77 - 115	203	183 - 223
IgA	176	141 - 211	427	384 - 470
IgG	855	684 - 1,026	2,149	1,934 -
2,364				
IgM	77	62 - 92	185	166 - 204
Prealbumin	20	16 - 24	41	37 - 45
Transferrin	238	190 - 286	455	409 - 501

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Тип: Control

KIT COMPOSITION:

Level 1	Human Plasma	5 x 0.5 mL, lyophilized
Level 2	Human Plasma	5 x 0.5 mL, lyophilized

INTENDED USE:

The Coagulation Control is intended to be used as a consistent test sample of known concentration for monitoring the performance of coagulation assays for the listed constituents.

CONTROL PREPARATION:

1. Reconstitute 1 vial of each Coagulation Control level with 0.5 mL of purified water and allow to sit for 30 minutes.
2. Gently swirl vial until dissolution is complete and the solution is homogenous (do not shake).
3. Replace the cap immediately after use to prevent evaporation or contamination.
4. Return the vial to 2-8°C promptly after use.

STORAGE AND HANDLING:

Store lyophilized and reconstituted controls at 2-8°C. Return all controls to 2-8°C promptly after use. Unopened controls can be used for 2 years from the date of manufacturing as indicated by the expiration date on the package and bottle labels. Reconstituted controls can be used for up to 7 days if stored at 2-8°C.

Reconstituted controls can be frozen one time and stored at -20°C for up to 5 weeks. Do not refreeze a second time.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Coagulation Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

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Intended Use

For in vitro diagnostic use only

The Control Set is to be used for monitoring the accuracy and precision of various CreatinKinase-MB assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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Intended Use

For *in vitro* diagnostic use only.

The Control Set is to be used for monitoring the accuracy and precision of various Lipids assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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The assay values and expected ranges are target values derived from inter laboratory data. The expected range values include variations of instrument and laboratory handling. The assay values were obtained using reagents available at the time of testing. Updates to the listed values may be made based upon additional data that becomes available or, if necessitated by a modification to a test method. The mean values established for your laboratory should fall within the ranges shown in Expected Results, however laboratory means may vary during the life of the control. Each laboratory should establish its own mean and precision parameters.

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Materials Provided

Control material.

Materials Required but not Provided

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Intended Use

For in vitro diagnostic use only.

The Control Set is to be used for monitoring the accuracy and precision of various HDL/LDL assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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Intended Use

For in vitro diagnostic use only.

The Control Set is to be used for monitoring the accuracy and precision of various Microalbumin assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

Product Description

The control material is supplied as a ready-to-use liquid requiring no reconstitution or dilution. It is prepared in an aqueous base fortified with reagent grade chemicals. Preservatives have been added to inhibit growth.

Precautions

Normal precautions exercised in handling laboratory reagents and/or infectious biological materials should be followed. Dispose of waste observing all local, state and federal laws.

Storage and Stability

Ensure that you tightly seal the vials after reconstitution and use to prevent evaporation during storage. Ensure that you store the vials upright to prevent spills or leakage.

	Storage	Stability
Unopened	2-8°C	Refer to the label for the expiration date.
Open Vial Stability	2-8°C	20 days

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Gently invert the vial and swirl to assure contents are homogenous.

Remove the cap from the vial.

Perform assay.

Record the results according to your quality assurance program.

Expected Results

Refer to the expected results supplied for mean and range value assignments. Verify vial lot number with those on the assay values sheet.

The assay values and expected ranges are target values derived from inter laboratory data. The expected range values include variations of instrument and laboratory handling. The assay values were obtained using reagents available at the time of testing. Updates to the listed values may be made based upon additional data that becomes available or, if necessitated by a modification to a test method. The mean values established for your laboratory should fall within the ranges shown in Expected Results; however laboratory means may vary during the life of the control. Each laboratory should establish its own mean and precision parameters.

Limitations

The results obtained using the quality control material depends on several factors: erroneous results can occur from improper storage, inadequate mixing, or sample handling errors associated with instrument or assay procedures. Do not use the quality control material if there is visible evidence of microbial growth in the vial. For more information about procedural limitations, refer to your instrument manual or assay product insert.

Disposing of Materials

Dispose of hazardous or biologically contaminated materials according to your institution's practices. Discard all materials in a safe and acceptable manner that is in compliance with all country, state, and local requirements.

Precautions

For in vitro diagnostic use only.

Handle serum and blood specimens as potentially infectious samples and follow the guidelines established by the Centers for Disease Control (CDC), Atlanta, Georgia, for blood collection and handling.

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Intended Use

For in vitro diagnostic use only

The Control Set is to be used for monitoring the accuracy and precision of various Microprotein assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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Storage	Stability
Unopened Open Vial Stability	Refer to the label for the expiration date. 20 days

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INTENDED USE:

The RF / ASO Liquid Control is intended for use as an assayed quality control material for serum protein analysis of Rheumatoid Factor and Anti-Streptolysin O.

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STORAGE AND HANDLING:

Store controls at 2-8°C. Return all controls to 2-8°C promptly after use. Unopened controls can be used up to the expiration date on the package and bottle labels.

Opened controls can be used for 60 days if stored tightly capped at 2-8°C.

EXPECTED VALUES:

The following values are given as an example. The expected values for the RF / ASO Liquid Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

ANALYTE MEAN (IU/mL)	ANALYZER RANGE (IU/mL)	LEVEL 1 MEAN (IU/mL)	LEVEL 2 RANGE (IU/mL)
ASO	Hitachi 911 103	82 - 124	295 236 - 355
RF	Hitachi 911 24	19 - 29	55 44 - 66
	Beckman CX7 27	22 - 33	72 57 - 86

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INTENDED USE:

The High-Sensitive CRP Control is intended for use as an assayed quality control material for monitoring the performance of C-Reactive Protein (CRP) immunoturbidimetric assays. FOR IN VITRO DIAGNOSTIC USE.

CONTROL PREPARATION:

1. Gently invert the bottle before use.
2. The control is ready to use.
3. Replace the cap immediately after use to prevent evaporation or contamination.

STORAGE AND HANDLING:

Store controls at 2-8°C. Return all controls to 2-8°C promptly after use. Controls can be used up to the expiration date shown on the vial label and package insert.

EXPECTED VALUES:

The following values are given as an example. The expected values for the High-Sensitive CRP Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

ANALYTE	ANALYZER	LEVEL 1	LEVEL 2
MEAN (mg/L)	RANGE (mg/L)	MEAN (mg/L)	RANGE (mg/L)
CRP (2) + Cal B	Hitachi 911 0.47	0.37 - 0.56	1.62
	Beckman CX7 0.55	0.44 - 0.66	1.62
1.30 - 1.95			
CRP (3) + Cal D	Hitachi 911 0.69	0.55 - 0.83	2.42
1.94 - 2.91			
CRP (3) + Cal E	Hitachi 911 0.54	0.43 - 0.64	1.72
1.38 - 2.07			
CRP (3) + Cal F	Hitachi 911 0.8	0.5 - 0.7	1.9
1.5 - 2.3			

Note: The above are standardized to the International Federation of Clinical Chemistry (IFCC) International Reference Preparation for Plasma Protein CRM 470.

Controls assayed with reagents from various manufacturers:

ANALYTE	ANALYZER	LEVEL 1	LEVEL 2
MEAN (mg/L)	RANGE (mg/L)	MEAN (mg/L)	RANGE
MEAN (mg/L)			
Hitachi 911	0.44	0.35 - 0.53	1.40
Beckman Coulter Image	0.53	0.42 - 0.64	2.11
1.69 - 2.53			
Olympus AU 400 / 640	0.53	0.42 - 0.64	1.65
1.32 - 1.98			

Precautions

For in vitro diagnostic use only. Handle serum and blood specimens as potentially infectious samples and follow the guidelines established by the Centers for Disease Control (CDC), Atlanta, Georgia, for blood collection and handling.

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Тип: Control

KIT COMPOSITION 1:

Level 1	Human Serum	1 x 5 mL, liquid (frozen)
Level 2	Human Serum	1 x 5 mL, liquid (frozen)
Level 3	Human Serum	1 x 5 mL, liquid (frozen)

KIT COMPOSITION 2:

Level 1	Human Serum	2 x 5 mL, liquid (frozen)
---------	-------------	---------------------------

KIT COMPOSITION 3:

Level 2	Human Serum	2 x 5 mL, liquid (frozen)
---------	-------------	---------------------------

KIT COMPOSITION 4:

Level 3	Human Serum	2 x 5 mL, liquid (frozen)
---------	-------------	---------------------------

INTENDED USE:

The Ferritin / IgE / Insulin Control is intended for use in the clinical laboratory as a consistent test sample of known concentration for monitoring assay conditions in many immunological determinations. Include Ferritin / IgE / Insulin Control with patient serum specimens when assaying for any of the listed constituents. Assay values are provided for the specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and instrument.

FOR IN VITRO DIAGNOSTIC USE

CONTROL PREPARATION:

1. Thaw at room temperature on a rocker or with periodic gentle inversion.
2. Once thawed, immediately store at 2-8°C. ONCE THAWED, DO NOT REFREEZE THE CONTROL.

STORAGE AND HANDLING:

Unopened vials of Ferritin / IgE / Insulin Control are stable for 90 days from receipt when stored at 2-8°C. Once opened, vials of Ferritin / IgE / Insulin Controls are stable for 30 days when stored tightly capped at 2-8°C. Ferritin / IgE / Insulin Control maintained at -20°C is stable until the expiration date on the box.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Ferritin / IgE / Insulin Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

ANALYTE	ANALYZER	UNITS	LEVEL 1 MEAN	LEVEL 2 RANGE	LEVEL 3 MEAN
Ferritin	Hitachi 911		17	14 - 21	153
	123 - 184		361	289 - 434	
	Beckman CX-7		19	16 - 23	132
	105 - 158		293	234 - 351	
Total IgE	Hitachi 911		31	25 - 38	145
	116 - 174		296	237 - 355	
Insulin	Hitachi 911		6.2	4.9 - 7.4	74.8
	59.9 - 89.8	104.2	83.4 - 125.0		
	Beckman CX-7		4.9	3.9 - 5.9	77.7
	62.2 - 93.3		>99	N/A	

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Biochemistry analyzer.

Analysis

Follow the QC procedure recommended by the instrument manufacturer.

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Тип: Control

KIT COMPOSITION:

Level 1 Human Plasma 5 x 1 mL, lyophilized
Level 2 Human Plasma 5 x 1 mL, lyophilized

INTENDED USE:

The Fibrinogen Control is intended for use in the quality control of fibrinogen assays.

FOR IN VITRO DIAGNOSTIC USE.

CONTROL PREPARATION:

1. Reconstitute with 1.0 mL of distilled water. Swirl gently and let stand undisturbed for 15 minutes at room temperature.
2. Do not invert vial or mix vigorously.

STORAGE AND HANDLING:

Reconstituted Level 2 Fibrinogen Control is stable for 8 hours when stored in capped vial at 2-8°C and Level 1 Fibrinogen Control is stable for 16 hours when stored in capped vial at 2-8°C. A precipitate may form when refrigerated. Gently warm the plasma to 37°C to minimize any precipitate.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Fibrinogen Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

METHOD	UNIT	LEVEL 1 RANGE	LEVEL 2 MEAN RANGE
Fibrinogen on Hitachi 911	MEAN		
548		146	117 - 175 457 365 -
Fibrinogen on Beckman Synchron CX	mg/dL	167	134 - 201 447 358 -
537			
Clauss Fibrinogen	mg/dL	117	N/A 524 N/A

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Intended Use

For in vitro diagnostic use only.

The Control Set is to be used for monitoring the accuracy and precision of various Fructosamine assay methods and to validate quantification of patient samples. The controls contain components of known concentrations and are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

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Expected Results

Refer to the expected results supplied for mean and range value assignments. Verify vial lot number with those on the assay values sheet. The assay values and expected ranges are target values derived from inter laboratory data. The expected range values include variations of instrument and laboratory handling. The assay values were obtained using reagents available at the time of testing. Updates to the listed values may be made based upon additional data that becomes available or, if necessitated by a modification to a test method. The mean values established for your laboratory should fall within the ranges shown in Expected Results; however laboratory means may vary during the life of the control. Each laboratory should establish its own mean and precision parameters.

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INTENDED USE:

The Cystatin C Control is intended for use as an assayed quality control material for monitoring the performance of cystatin C assays.
FOR IN VITRO DIAGNOSTIC USE.

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STORAGE AND HANDLING:

Store controls at 2-8°C. Return all controls to 2-8°C promptly after use. Unopened controls can be used for one year from the date of manufacturing as indicated by the expiration date on the package and bottle labels.

Opened controls can be used for 1 month if stored at 2-8°C.

EXPECTED VALUES:

The following values are given as an example. The expected values for the Cystatin C Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert and/or vial labels for each lot for the expected values.

REAGENT	LEVEL 1	LEVEL 2
MEAN (mg/L)	RANGE (mg/L)	MEAN (mg/L) RANGE (mg/L)
Cystatin C	0.8 0.7 - 0.9	1.6 1.4 - 1.8

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Intended Use

For *in vitro* diagnostic use only.
The Chemistry Control is to be used for monitoring the accuracy and precision of clinical chemistry procedures. This control contains constituents commonly of interest in a general chemistry control, including common drugs and thyroid. The assayed multi-analyte control product is packaged with one level of analytes.

Product Description

The quality control material is prepared from human serum with enzymes, nonprotein constituents, non-human protein, and bacteriostatic agents added. The constituents were adjusted to the levels listed in Expected Values.

Precautions

BIOHAZARD: Human source material. Handle as if potentially infectious. Human serum was used in the manufacture of this product. Each donor unit used was tested by FDA-approved methods and found nonreactive for hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), and HIV 1 and HIV 2. Because no known test method can offer complete assurance that infectious agents are absent, all products containing human source material should be handled in accordance with recommendations from Centers for Disease Control/National Institute of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories, 1999." This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with large volumes of water to prevent azide build-up.

Storage and Stability

Ensure that you tightly seal the vials after reconstitution and use to prevent evaporation during storage. Ensure that you store the vials upright to prevent spills or leakage. CK and Bilirubin are sensitive to light. Store the vials away from light.

Storage Stability

Unreconstituted 2-8°C Refer to the label on each vial and on the package for the expiration date.

Reconstituted 2-8°C 7 days

Exceptions: Bilirubin and Alkaline Phosphatase which are stable 48hrs. ALP may increase with time.

General Instructions for Use

Use the quality control material according to the directions accompanying the instrument or the assay procedure used. Treat the quality control material in the same manner as patient samples.

1. Remove the screw cap and gently remove the rubber stopper from the vial.
2. Pipette exactly 5.0 mL of distilled or deionized water to the vial using a volumetric pipette.
3. Replace the stopper in the vial, allow the vial to sit for 10 minutes.
4. Gently invert the vial three (3) times and swirl until the contents are homogenous.
5. Record the results according to your quality assurance program.

Expected Results

Refer to the Expected Values table supplied for assay mean and ranges. Before use, verify that the vial lot number corresponds to the lot number listed on the Expected Values table.

The Expected Values and ranges are target values derived from inter laboratory data. The expected range values include variations of instrument and laboratory handling. The assay values were obtained using in-date HTI reagents available at the time of testing. Updates to the listed values may be made based upon additional data that becomes available or, if necessitated by a modification to a test method. The mean values established for your laboratory should fall within the ranges shown in Expected Values, however laboratory means may vary during the life of the control. Each laboratory should establish its own mean and precision parameters.

Limitations

The results obtained using the quality control material are dependant upon several factors: erroneous results can occur from improper storage, reconstitution errors, inadequate mixing, or sample handling errors associated with instrument or assay procedures. Do not use the quality control material if there is visible evidence of microbial growth in the vial or a lack of vacuum when opening the vial for the first time. For more information about procedural limitations, refer to your instrument manual or assay product insert.

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Intended Use

The product is for use as a calibrator of HTI Inc. clinical chemistry assays. This calibrator material is well suited for automated and semi-automated analytical procedures.

Summary

The HTI, Inc. calibrator is a human based serum. The concentration of the calibrator components have been adjusted to ensure optimal calibration of the HTI, Inc. methods on specified analyzers.

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Precautions

For *in vitro* diagnostic use only: Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2, 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens. (29CFR part 1910.1030). Federal register, July 1, 1998, 6:287-289.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For *in vitro* diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2, 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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1. Accurate volumetric pipetting devices.
2. Timer.
3. General laboratory equipment.

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998. 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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High Technology, Inc

Intended Use

The product is for use as a calibrator of HTI inc. clinical chemistry assays. This calibrator material is well suited for automated and semi-automated analytical procedures.

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Handling Instructions

Carefully open one bottle, avoiding the loss of lyophilizate. Using a volumetric pipette, add exactly 5.0 ml of diluent to the lyophilized serum. Gently invert the vial intermittently over a period of 20 minutes to ensure complete dissolution of contents. Immediately prior to use, gently invert the vial 5 -10 times.

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Summary

The HTI, Inc. calibrator is a human based serum.

The concentration of the calibrator components have been adjusted to ensure optimal calibration of the HTI, Inc. methods on specified analyzers.

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Calibrator Values¹

The calibrator values were determined using HTI, Inc. reagent methods and the analyzers listed in the value assignment table (See reverse side of package insert.) Determinations were performed under strictly standardized conditions, utilizing known reference materials. Traceability information available upon request.

Calibrator Storage and Stability

Unreconstituted chemistry calibrator is stable until the expiration date when stored at 2-8°C. Reconstituted chemistry calibrator is stable for seven days when stored at 2-8°C with the exception of Bilirubin, which is stable five days at 2-6°C. Store calibrator tightly capped and protected from light when not in use.

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References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030). Federal register July 1, 1991. 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For *in vitro* diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2. 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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Assay

Follow the calibration procedure recommended by the instrument manufacturer

References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998. 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from Oct. 17th, 2000.

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Precautions

For in vitro diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV-1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2. Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998, 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For in vitro diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HCV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2. Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030) Federal register. July 1, 1998. 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Calibrator Values

The calibrator values were determined using HTI Inc. reagent methods and the analyzers listed in the value assignment table. (See reverse side of package insert.) Determinations were performed under strictly standardized conditions, utilizing known reference materials. Traceability information available upon request.

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998. 6:267-280
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Handling Instructions

Carefully open one bottle, avoiding the loss of lyophilizate. Using a volumetric pipette, add exactly 5.0 ml of diluent to the lyophilized serum. Gently invert the vial intermittently over a period of 20 minutes to ensure complete dissolution of contents. Immediately prior to use, gently invert the vial 5 - 10 times.

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References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998, 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For *in vitro* diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV, and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2, 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998. 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For *in vitro* diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV-1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2. 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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1. Accurate volumetric pipetting devices
2. Timer
3. General laboratory equipment

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens (29CFR part 1910.1030). Federal register July 1, 1998, 6:267-280
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Summary

The HTI Inc. calibrator is a human based serum. The concentration of the calibrator components have been adjusted to ensure optimal calibration of the HTI Inc. methods on specified analyzers.

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Precautions

For *in vitro* diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV-1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2, 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030). Federal register. July 1, 1998, 6:287-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For in vitro diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens (29CFR part 1910.1030). Federal register. July 1, 1998; 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For in vitro diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2. 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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Summary

The HTI Inc. calibrator is a human based serum. The concentration of the calibrator components have been adjusted to ensure optimal calibration of the HTI, Inc. methods on specified analyzers.

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Calibrator Storage and Stability

Unreconstituted chemistry calibrator is stable until the expiration date when stored at 2-8°C. Reconstituted chemistry calibrator is stable for seven days when stored at 2-8°C with the exception of Biliubin, which is stable five days at 2-8°C. Store calibrator tightly capped and protected from light when not in use.

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References

1. Department of Labor, Occupational Safety and Health Standards. Bloodborne pathogens. (29CFR part 1910.1030) Federal register. July 1, 1998, 6:267-280
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Precautions

For in vitro diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV-1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2, 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

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References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens. (29CFR part 1910.1030). Federal register, July 1, 1998, 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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Intended Use

The Multi-Analyte Chemistry calibrator is for use as a calibrator of HTI Inc. clinical chemistry assays. This calibrator material is well suited for automated and semi-automated analytical procedures.

Summary

The HTI, Inc. Multi-Analyte calibrator is a human based serum. The concentration of the calibrator components have been adjusted to ensure optimal calibration of the HTI, Inc. methods on specified analyzers.

Product Description

The product consists of lyophilized human serum and a diluent for reconstitution. The human serum contains additives to provide the defined assay values. The concentrations of the calibrator components are lot specific.

Calibrator Values

The calibrator values were determined using HTI, Inc. reagent methods and the analyzers listed in the value assignment table. (See reverse side of package insert.) Determinations were performed under strictly standardized conditions, utilizing known reference materials. Traceability information available upon request.

Calibrator Storage and Stability

Unreconstituted chemistry calibrator is stable until the expiration date when stored at 2-8°C. Reconstituted chemistry calibrator is stable for seven days when stored at 2-8°C with the exception of Bilirubin, which is stable five days at 2-8°C. Store calibrator tightly capped and protected from light when not in use.

Precautions

For in vitro diagnostic use only. Human serum was used in the manufacture of this product. Each donor unit was tested for antibodies to HIV1/2, HCV and found to be non-reactive for HBsAg and HIV-1Ag by FDA accepted test methods. Because no test method can offer complete assurance that products derived from blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens. In the event of exposure, the directives of the responsible health authorities should be followed. 2, 3 Safety data sheets are available upon request. Disposal of all waste material should be in accordance with local guidelines.

Handling Instructions

Carefully open one bottle, avoiding the loss of lyophilizate. Using a volumetric pipette, add exactly 5.0 ml of diluent to the lyophilized serum. Gently invert the vial intermittently over a period of 20 minutes to ensure complete dissolution of contents. Immediately prior to use, gently invert the vial 5-10 times.

Materials Provided

Multi-Analyte calibrator with Diluent

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1. Accurate volumetric pipetting devices
2. Timer
3. General laboratory equipment.

Assay

Follow the calibration procedure recommended by the instrument manufacturer.

References

1. Department of Labor, Occupational Safety and Health Standards: Bloodborne pathogens. (29CFR part 1910.1030). Federal register: July 1, 1998. 6:267-280.
2. Council Directive (2000/54/EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.

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STANDARD

For monitoring the accuracy and precision of assay procedures used commonly in clinical chemistry.
For *in vitro* diagnostic use.

INTRODUCTION

Product contains a broad range of chemicals in a general chemistry calibrator. The human serum calibrator is to be used in the same manner that an unknown specimen would be used in accordance with the directions of the assay procedure.

REAGENT COMPOSITION

Product is prepared from human serum with human and nonhuman enzymes and non-protein constituents added.

WARNINGS AND PRECAUTIONS

For *In-Vitro* Diagnostic Use Only.

Human source material was used in the manufacture of this product. Each donor unit used was tested with licensed reagents and found negative for antibodies to HIV1/2 and HCV and non-reactive for HBsAg and HIV-1 Ag by FDA-accepted methods. Because no test method can offer complete assurance that HIV, hepatitis B virus, HCV or other infectious agents are absent, this material should be handled as if capable of transmitting infectious diseases.

REAGENT STORAGE

1. Reagent is stable until the expiration date on the label when stored tightly capped in amber glass vial at 2 - 8°C.
2. Product is stable for 5 days after open when stored at 2 - 8°C.
3. Do not use after the expiration date.

REAGENT DETERIORATION

1. Failure to obtain accurate results in the assay of control materials may indicate reagent deterioration.
2. If there is visible evidence of microbial growth or excessive turbidity in the product, do not use that vial. Discard the vial according to the requirements of your local waste management authorities.

MATERIALS PROVIDED

Standard

EXPECTED VALUES

Listed values were obtained using multiple lots of reagents, in date at the time of testing, with various analyzers. End-user is responsible for determining if values are appropriate to their system in use.

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RESULTS REGISTRATION

See User Manual for your Biochemistry analyzer.

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WARNINGS AND PRECAUTIONS

For *In-Vitro* Diagnostic Use Only.

Human source material was used in the manufacture of this product. Each donor unit used was tested with licensed reagents and found negative for antibodies to HIV-1/2 and HCV and non-reactive for HBsAg and HIV-1 Ag by FDA-accepted methods. Because no test method can offer complete assurance that HIV, hepatitis B virus, HCV or other infectious agents are absent, this material should be handled as if capable of transmitting infectious diseases.

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High Technology, Inc

STANDARD

For monitoring the accuracy and precision of assay procedures used commonly in clinical chemistry. For *in vitro* diagnostic use.

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REAGENT COMPOSITION

Reagent is prepared from human serum with human and nonhuman enzymes and non-protein constituents added.

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REAGENT STORAGE

1. Product is stable until the expiration date on the label when stored tightly capped in amber glass vial at 2 - 8°C.
2. Product is stable for 5 days after open when stored at 2 - 8°C.
3. Do not use after the expiration date.

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