Siemens Healthcare GmbH, SHS AT QT RA, Siemensstr. 1, 91301 Forchheim

To the Federal service for surveillance in
Healthcare in the Russian Federation

Name
Department

Karla Flesner-Oeste
SHS AT QT RA

Mobile
E-mail

+49 (1525) 9518360
karla.flesner-oeste@siemens-healthineers.com

Date

January 24, 2022

Declaration Letter Operator Manual

We, Siemens Healthcare GmbH, hereby confirm that the following Operator Manual for Cios Fusion has been printed by the manufacturer in Germany:

XPR2-370G.620.01.02.02

Sincerely yours,

Siemens Healthcare GmbH

Siegfried Quinger
Head of Quality
Business Area Advanced Therapies

Ali-Ihsan Can
Head of Regulatory Affairs
Business Area Advanced Therapies

Siemens Healthcare GmbH
Management: Bernhard Montag, President and Chief Executive Officer;
Darleen Caron, Jochen Schmitz, Christoph Zindel

Chairman of the Supervisory Board: Ralf P. Thomas
Registered office: Munich, Germany; Commercial Registry: Munich, HRB 213821
WEEE-Reg.-No. DE 64872105

Siemensstr. 1
91301 Forchheim
Germany

Tel.: +49 (9191) 18 0
Fax: +49 (9191) 18 9999
siemens.com/healthcare

SIEMENS

www.siemens.com/healthcare

Cios Fusion

Operator Manual

Cios Fusion

Operator Manual

Legend

	Indicates a hint Is used to provide information on how to avoid operating errors or information emphasizing important details
	Indicates the solution of a problem Is used to provide troubleshooting information or answers to frequently asked questions
	Indicates a list item
	Indicates a prerequisite Is used for a condition that has to be fulfilled before starting a particular operation
	Indicates a one-step operation
	Indicates steps within operating sequences
<i>Italic</i>	Is used for references and for table or figure titles
	Is used to identify a link to related information as well as previous or next steps
Bold	Is used to identify window titles, menu items, function names, buttons, and keys, for example, the Save button
Blue	Is used to emphasize particularly important sections of the text
UI text	Is used for on-screen output of the system including code-related elements or commands
User text	Is used to identify inputs you need to provide
Menu > Menu Item	Is used for the navigation to a certain submenu entry
 CAUTION	CAUTION Used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury or material damage. CAUTION consists of the following elements: ■ Information about the nature of a hazardous situation ■ Consequences of not avoiding a hazardous situation ■ Methods of avoiding a hazardous situation
 WARNING	WARNING Indicates a hazardous situation which, if not avoided, could result in death or serious injury. WARNING consists of the following elements: ■ Information about the nature of a hazardous situation ■ Consequences of not avoiding a hazardous situation ■ Methods of avoiding a hazardous situation

1	Introduction	7
1.1	Use of the Product	8
1.2	Information about this Operator Manual	10
2	Safety	11
2.1	General Safety Information	13
2.2	Personal Safety	23
2.3	Equipment Safety	42
3	System Description	49
3.1	Device Description	51
3.2	System Operation	63
3.3	Emergency Situations	90
4	Examination	95
4.1	Safety information	97
4.2	Preparing the system	99
4.3	Patient registration	103
4.4	Defining the examination settings	114
4.5	Acquiring images	123
4.6	Ending the examination	131
4.7	Performing special examinations	133
5	Postprocessing	141
5.1	Patient data	143
5.2	2D image processing	152
5.3	Subtraction processing	176
6	Documentation	185
6.1	Filming/Printing	186
6.2	Exporting	189
6.3	Checking the data transfer	199
6.4	Reports	201
7	Configuration	205
7.1	User Settings	207
7.2	Examination settings (PEX Editor)	219
8	Maintenance	233
8.1	Functional and safety checks	234
8.2	Service via network connection	239
8.3	Cleaning and Disinfection	244
8.4	Protection classes/Protection measures	246
8.5	Active ingredient classes	246
9	Technical Specifications	247
9.1	Curves and diagrams	249
9.2	Notes on electromagnetic compatibility (EMC)	261

Table of contents

9.3	Original equipment manufacturer	265
9.4	System data	266
9.5	List of Standards	275
9.6	Labels	276
■	10 Options	287
10.1	System accessories	288

1 Introduction

1.1	Use of the Product	8
1.1.1	Intended use	8
	Body region	8
1.1.2	Patient group	8
1.1.3	Contraindications	8
1.1.4	Physical functionality	8
1.1.5	User profile	8
	Training	9
1.1.6	Conditions of use	9
	Maintenance, cleaning and disinfection, service	9
1.1.7	Essential performance characteristics	9
1.1.8	Frequently-used operating functions	9
1.1.9	Operating functions regarding safety	9
1.2	Information about this Operator Manual	10
1.2.1	Names and parameters	10
1.2.2	Trademarks	10
1.2.3	Values	10
1.2.4	Layout conventions	10

1.1 Use of the Product

We welcome you as a user of the powerful Cios Fusion, the multifunctional C-arm system from Siemens.

1.1.1 Intended use

The Cios Fusion mobile radiography system was developed for clinical applications such as vascular, cardiac, traumatological, orthopedic, neurological, endoscopic, urologic, and pain therapy procedures, as well as cholangiography and pediatric examinations. It supports medical examinations through digital imaging.

Body region

The systems are suitable for a large variety of clinical visualizations such as

- Visualization of bones
- Visualization of vessels, soft tissues
- Visualization of implants
- Visualization of surgical instruments

1.1.2 Patient group

In principle, the Cios Fusion can be used for examining all types of patients without restriction. When using this system, the relevant applicable country-specific requirements must be complied with.

1.1.3 Contraindications

Prior to the examination, it must be confirmed by a physician that an examination is permissible and checked whether increased precautions are necessary.

1.1.4 Physical functionality

The Cios Fusion X-ray system is a compact and powerful mobile X-ray system with flat detector and digital image processing.

The system is equipped with a footswitch and a hand switch for the release of radiation.

The following modes of operation are available in Cios Fusion for the broad scope of applications: Single image radiography (Single Image), fluoroscopy (Fluoro), subtraction (Sub), and roadmap (Road).

1.1.5 User profile

The system can be used properly only by suitably qualified medical personnel. This requires operating personnel to be familiar with the Operator Manual. This manual must be studied in detail prior to starting up the system. Special attention must be given to the following sections:

- General Safety Information
- Personal safety

- Equipment safety
- Maintenance

It is the responsibility of the system operator to ensure that operating personnel are professionally and properly instructed.

As the manufacturer of the system, SIEMENS is not responsible for any impact on its safety, reliability and performance if it is not used in accordance with the Operator Manual.

Please also read the supplements and addenda to the Operator Manual necessitated by technical developments.

Training

It is recommended that operating personnel undergo application training prior to using the system. This training can be conducted by Siemens employees.

1.1.6 Conditions of use

The digital Cios Fusion X-ray system is a mobile system that is designed for use in the ER, OR, and in radiology and endoscopy departments of hospitals, clinics, and outpatient practices. The applicable country-specific requirements must be observed when using the system.

The system is suited for operation on hospital networks (IEC 60601-1 Chap. 14.13 IT security statement). For more information, see the technical documentation.

Maintenance, cleaning and disinfection, service

See (→ Page 233 *Maintenance*).

1.1.7 Essential performance characteristics

See (→ Page 247 *Technical Specifications*).

1.1.8 Frequently-used operating functions

- System On/Off
- Patient registration
- Selecting the operating mode
- Acquiring X-ray images
- Processing X-ray images (Postprocessing/Adding comments)
- Saving/archiving studies locally or to a PACS

1.1.9 Operating functions regarding safety

- Emergency stop function (Emergency stop button)
- Radiation protection (Collimator function)
- Moving the C-arm

1.2 Information about this Operator Manual

1.2.1 Names and parameters

All names and data on patients and equipment that are used as examples in this Operator Manual are entirely fictional.

Any resemblance to names of real persons and institutions is entirely coincidental.

All parameters and images shown in this Operator Manual are examples. Only the parameters displayed by your system are definite.

1.2.2 Trademarks

Transliner® is a registered trademark of Siemens AG.

Microsoft and Adobe Acrobat are registered trademarks, and Windows is a trademark of Microsoft Corporation.

Oracle and Java are trademarks or registered trademarks of Oracle America, Inc.

All other product or company names mentioned in this document are trademarks or registered trademarks of their respective owners and are used only for purposes of identification or description.

1.2.3 Values

All numbers specified are typical values unless specific tolerances are indicated.

1.2.4 Layout conventions

The Operator Manual has several registers. A detailed table of contents listing all chapters contained in the manual is provided at the beginning.

Certain sections of text are marked with symbols to help you quickly identify the information content of the text.

See (→ Page 4 *Legend*).

2 Safety

2.1	General Safety Information	13
2.1.1	Laws and regulations	13
2.1.2	Range of application	13
2.1.3	Using the system	13
2.1.4	Software	14
2.1.5	Potential equalization	14
2.1.6	Electromagnetic compatibility	15
2.1.7	Use in connection with high frequency	15
2.1.8	Maintenance and inspection	15
	Wear and tear	15
	Image quality	16
	Calibration	16
	Performing maintenance	16
2.1.9	Malfunctions	16
	Error messages at the C-arm system	16
	System messages on the monitor	17
	Handling error messages	18
2.1.10	Malfunction of electrical systems	18
	System failure	19
	Switching to emergency power supply	19
	Disconnecting the power plug	20
2.1.11	EMERGENCY STOP	20
	Press the EMERGENCY STOP button	20
	Unlocking the EMERGENCY STOP button	21
2.1.12	Fire protection	21
2.1.13	Explosion protection	22
2.1.14	Overload protection	22
2.2	Personal Safety	23
2.2.1	Open heart and skull examinations	23
2.2.2	Crushing hazards on the C-arm system	23
2.2.3	Mechanical damage	24
2.2.4	Radiation protection	24
	For the patient	25
	For the operating personnel	25
	For patients and operating personnel	25
	Deterministic radiation effect	26
	Location and size of the relevant operating areas	27
	Radiation interruption for all operating modes	41
2.3	Equipment Safety	42
2.3.1	Mechanical safety	42
2.3.2	Positioning the C-arm	42
	Brakes	42
	Transport	42
2.3.3	Installation, repair	43
2.3.4	Original accessories	43
2.3.5	Combination with other systems	44

2 Safety

2.3.6	Attachment of dedicated options	44
	General safety requirements	44
	Attenuation equivalent	44
	Weight counterbalance	45
	Image quality	45
	Electrical safety	45
	Electromagnetic compatibility	45
	Attaching navigation systems	45
	Additional safety information	47
2.3.7	Disposal	47

2.1 General Safety Information

2.1.1 Laws and regulations

If legally binding regulations govern the installation and/or operation of the system, it is the responsibility of the installer and/or the operator to observe these regulations.

Regulations required by law and the radiation protection regulations must be observed in all countries. Deviating from this Operator Manual, values may be set according to country-specific regulations.

The product bears a CE mark in accordance with the provisions of Directive 93/42/EEC dated June 14, 1993 on medical devices as well as Directive 2011/65/EU dated June 8 2011 on the restricted use of certain hazardous substances in electrical and electronic devices.

Data related to individual persons are subject to data protection. Please comply with the applicable legal regulations.

Legally required tests must be performed at the specified intervals. These tests include, for example,

- Constancy test according to the X-ray ordinance (§16 RöV) in the Federal Republic of Germany.
- Tests based on DHHS guidelines (Department of Health and Human Services) where applicable.

2.1.2 Range of application

This Operator Manual is valid for the following product:

- Cios Fusion

2.1.3 Using the system

CAUTION

Improper use of the system beyond its intended use.

Risk to the patient and user!

- ◆ The system may only be used as defined by its intended use for clinical applications.

CAUTION

System operation by untrained users.

Risk of improper operation due to misinterpretation of image information.

- ◆ Only persons with the required knowledge and expertise who have undergone the appropriate user training may work with the system.

CAUTION

Occurrence of audible signals.

Hazards possible, for example due to crushing risks, extended radiation release, and high dose rate.

- ◆ Please note that continued system operation may lead to one of these hazards.
- ◆ Use the system with extreme caution to avoid these hazards.

2.1.4 Software

The system and user software used in this product is protected by copyright. The current software version is displayed on screen during system start-up.

CAUTION

Use of unreleased or altered software or hardware components.

Risk of malfunction that could endanger the patient or product.

- ◆ Only use Siemens-authorized software and hardware components.
- ◆ Only perform repairs, or have them performed by an authorized third party, after having received our express written approval.

CAUTION

Unauthorized or incorrect alterations/changes of the software or connection of the system to a network.

Risk of unauthorized access.

- ◆ Make sure that all necessary precautions are taken with regard to the applicable security level if a functionality is added or the configuration of the condition on delivery is changed.

CAUTION

Unauthorized changes or interventions in the imaging system.

No or improper function of the imaging system.

- ◆ Do not perform any changes or interventions of any type on the imaging system without the written consent of Siemens Healthcare (supplier). This includes replacing hardware or installing and running additional software.

2.1.5 Potential equalization

Products for which equipotential bonding is required may only be operated in medical facilities where supplementary equipotential bonding has been installed and tested according to DIN VDE 0107/10.94 section 4 for Germany or the relevant country-specific regulations.

2.1.6 Electromagnetic compatibility

This medical device complies with the requirements of the applicable standard on electromagnetic compatibility (EMC). See (→ Page 261 *Notes on electromagnetic compatibility (EMC)*)

Please be advised that other mobile electronic devices, e.g. cellular telephones, exceeding the established emissions limits in the EMC standard may disrupt the functions of your medical device.

CAUTION

Interference.

Impact on the patient's electronic life-sustaining systems.

- ◆ In case of unusual system behavior (performance characteristics), additional measures (adjustment, repositioning) may be required.

2.1.7 Use in connection with high frequency

The following regulations for use must be observed:

■ IEC/TR 1289-1/07.94/

High frequency surgical equipment - Part 1: Operation

■ IEC/TR 1289-2/08.94/

High frequency surgical equipment - Part 2: Maintenance

2.1.8 Maintenance and inspection

Before using the equipment for examination, the user must ascertain that all safety-relevant devices function properly and that the system is ready for operation.

Wear and tear

The system is subject to mechanical and electrical wear and tear. In the interest of the safety of patients, operating personnel and third persons, maintenance and safety checks must be carried out every 24 months to maintain the operational safety and reliability of the product.

CAUTION

Wear and tear.

**Risk of injury to the patient, operating personnel, and other persons.
Permanent damage to the system.**

- ◆ Follow the maintenance guidelines to maintain the safety and functionality of the system.

Please observe the relevant information in (→ Page 233 *Maintenance*).

Image quality

Maintenance should include checking the image quality. Maintenance at regular intervals is recommended to always ensure best image quality.

To ensure optimal image quality, have the following functions checked in particular as part of regular maintenance:

Pixel shift, image rotation, noise reduction, edge enhancement, subtraction, roadmapping, noise filter.

Calibration

Maintenance should include checking the flat detector calibration.

CAUTION

Geometric calibration adjusted incorrectly.

Risk of incorrect evaluation.

- ◆ Perform calibrations on a regular basis.

Performing maintenance

Maintenance work should be performed by trained technical personnel only. If you do not have a maintenance contract, please contact Siemens Customer Service.


If national laws or regulations specify more frequent inspection and/or maintenance, these must be observed.

2.1.9 Malfunctions

In the event of malfunctions of the Cios Fusion system, please call Siemens Customer Service.

Error messages at the C-arm system

When a malfunction is detected, the Cios Fusion system is disabled. The problem is shown as an error on all control units.


A message window with additional explanations is also displayed on the left monitor (→ Page 17 *System messages on the monitor*).

System messages on the monitor

On the left monitor (upper left), three different types of system messages can appear. The type of message is identified by a corresponding symbol (top).

Error message:

Warning:

Information:

Handling error messages

- ✓ An error message window is displayed.

Acknowledging error messages

- 1 Please read the error messages carefully. Use the scroll bars to display error messages outside the area currently displayed.
- 2 Acknowledge the error message by clicking OK or pressing the footswitch or hand switch.

If you cannot continue with the normal examination mode, contact Siemens Customer Service.

Repeatedly occurring errors

If errors occur repeatedly, switch off the Cios Fusion and notify Customer Service. Save the log file with the logged system activities beforehand:

- 1 Right-click this icon on the lower left of the left monitor.
- 2 In the dialog box, select the storage destination (USB and/or local storage for analysis via remote diagnostics) and confirm with **OK**.
- 3 Additionally, write down the following information:
 - Error number and time when error occurred.
 - Operating mode selected.
 - Was radiation activated when the error occurred?
 - Is the error related to an operating process?
- 4 In the case of a malfunction or failure of the radiation indicator, please notify Siemens Customer Service.

2.1.10 Malfunction of electrical systems

In case of risks for patients and operators (e.g. if there is no live image on the monitor and the radiation indicator is on despite this) or in case of risks for the device, you must disconnect the power plug immediately. Cios Fusion is completely shut down and disconnected from line power. This

- switches off radiation
- aborts the current system program
- aborts and cancels current operating sequences
- deletes all image information not saved to a hard disk

CAUTION

Technical error.

Unintentional radiation release.

- ◆ If the unintentional radiation release cannot be ended with the EMERGENCY STOP button, press the off switch to shut down the entire system.

Only after the cause of the hazard has been clearly identified and remedied may the system be reconnected to the power supply. In all other cases, e.g. system malfunction, contact Siemens Customer Service immediately.

CAUTION

Failure of LED indicators on the control panel.

Treatment not possible.

- ◆ Perform the daily function check of the LED indicators before beginning treatment.
- ◆ If an LED indicator fails, contact Siemens Customer Service.

System failure

The user must have a replacement unit available if a system failure could predictably cause a critical situation resulting in patient injury during a medical examination.

CAUTION

System component failure during cardiac/vascular procedure.

Risk of injury to the patient!

- ◆ Prepare a reserve system during cardiac/vascular applications.

CAUTION

Due to the complexity of the system or if the line voltage is outside of the designated specification, the loss of X-ray imaging or other system functions during an examination cannot be entirely excluded.

Imaging is delayed or not possible at all, resulting in suboptimal treatment.

- ◆ For this reason, establishing emergency procedures for such cases should be considered.

Switching to emergency power supply

In the event of a power interruption of longer than 10 ms, the Cios Fusion switches off. In this case the Cios Fusion must be switched on again after the system has switched over to the emergency power supply.

Disconnecting the power plug

After the power plug is disconnected, the uninterruptible power supply (UPS) supplies voltage to the imaging system and the left-hand monitor until the Cios Fusion switches off completely.

CAUTION

Unsuitable devices connected to the UPS.

Risk of electrical shock or damage to the system.

- ◆ Do not attach additional devices to the UPS sockets.

2.1.11 EMERGENCY STOP

Press the EMERGENCY STOP button

Immediately press the red EMERGENCY STOP button (arrow) on the C-arm system control unit at the first sign of a dangerous situation resulting from motorized movement.

- Motorized vertical movement is then immediately blocked.
- Active radiation is terminated.
- All other system functions remain unaffected by this.

CAUTION

Manual or motorized movements.

Risk of crushing and collision for persons and objects in the vicinity of moving parts.


- ◆ Check the EMERGENCY STOP function a daily basis.

CAUTION

Position of EMERGENCY STOP button not known.

In an emergency, the EMERGENCY STOP button will not be actuated in a timely manner.

- ◆ Make sure you know where the EMERGENCY STOP buttons are located. EMERGENCY STOP buttons are in the following locations (see the figure).

Emergency Stop button


If during an examination an error occurs and radiation cannot be interrupted by releasing the trigger switch, press the EMERGENCY STOP button on the chassis of the C-arm system.

Unlocking the EMERGENCY STOP button

Unlock the button only after the danger has clearly been eliminated.

- ◆ To unlock the button, gently turn it clockwise so that it pops back out.

2.1.12 Fire protection

CAUTION

Fire in or near the system.

Injury to patient and personnel, and damage to device. Risk of gas poisoning due to burning plastic.

- ◆ In case of fire switch off the system.
- ◆ Be aware of and inform patients of escape routes.
- ◆ Be aware of where fire extinguishers are located and know how to use them.

Please inform our Customer Service prior to starting up the Cios Fusion again if repair work has to be performed due to the fire.

2.1.13 Explosion protection

CAUTION

The system is not designed for operation in areas where there is an explosion hazard. The system does not meet the requirements for AP/APG classification.

Explosion hazard!

- ◆ The system may not be used in locations where explosive atmospheres may occur.

2.1.14 Overload protection

Prolonged continuous radiation at maximum tube load is possible in fluoroscopy mode. However, this can cause the X-ray tube assembly to heat up. For this reason, the X-ray tube assembly and C-arm have a thermal monitor. If the X-ray tube assembly is too warm, there may be a reduction in pulse rate if necessary in all operating modes beginning with the next scene. A corresponding message will be displayed.

CAUTION

Increased temperature during intensive use.

Risk of burns!

- ◆ The allowed temperature for brief contact is exceeded when the temperature switch of the single tank is triggered.

CAUTION

Defect in the radiation generation cooling system.

Risk of burns!

- ◆ Before beginning the examination, perform the daily function and safety checks.

CAUTION

Under extreme operating conditions (continuous operation), the X-ray tube can become very hot.

Risk of burns!

- ◆ Avoid touching the X-ray tube housing.

A constant X-ray tube load can lead to temperature-related cutoff of the radiation. Please make sure that the system has cooled down before performing critical interventional applications. We also recommend keeping a second system ready for continuous radiation lasting longer than 30 minutes.

In the rare event image display is delayed due to high system utilization, radiation is automatically turned off. To continue treatment, radiation release has to be actuated again.

If the Radiography mode is misused on purpose by the operator for real-time imaging, the image display delay may be longer than in Radioscopy.

2.2 Personal Safety

2.2.1 Open heart and skull examinations

If an approved system is used alone or with other equipment for open heart or open skull examinations, a conductive connection must be made between the system and a potential equalization point, e.g. the tabletop.

(→ Page 65 *Establishing the equipotential bonding connection*)

Only then can the patient be connected to the system.

2.2.2 Crushing hazards on the C-arm system

Correct handling of the C-arm system requires that operating personnel and patients use only the handles provided for this purpose. Where this is not possible, monitor the points of potential crush injury between movable system parts and their guide openings.

CAUTION

The distance between the patient or other persons and the system is too small.

Risk of crushing!

- ◆ Use the system only if the patient and other persons can be observed during movements.

CAUTION

The points marked on the figure indicate hazardous locations around the system.

Risk of injury to the patient and personnel due to crushing or collision.

- ◆ Be careful around the hazardous locations indicated.

CAUTION

The distance between the handles and other components is too small.

Risk of crushing!


- ◆ Be careful around the hazardous locations indicated.

CAUTION

Movements of the C-arm during patient positioning.

Crushing hazard for the patient!

- ◆ Be careful around the hazardous locations indicated.

Hazardous locations during C-arm movements

2.2.3 Mechanical damage

To avoid injury to the patient, operating personnel or third parties, mechanical damage to the system must be repaired by authorized service personnel.

2.2.4 Radiation protection

The unit is intended for procedures that may involve high skin doses when it is used as intended (mainly due to the long examination times), which presents a risk of radiation injury.

Therefore observe the following important notes in order to keep the dose absorbed by the patient as low as possible.

For the patient

- Keep the radiation field as small as possible without reducing the active measuring field.
- If possible, ensure the best possible protection of the patient during fluoroscopy and acquisitions in the vicinity of his or her reproductive organs (use gonadal shields, ovarian shields and lead rubber covers).
- Remove all radiopaque parts from the fluoroscopy field or field of view, if possible.
- Set the voltage of the X-ray tube as high as necessary (note the image quality, however).

For the operating personnel

- Wear protective clothing in the control area during an examination: Depending on the application involved, the use of lead vests, lead jackets, X-ray protection goggles, sternum protection and protective gloves is recommended.

For patients and operating personnel

- Keep the fluoroscopic time as short as possible.
- Set the tube-skin distance as high as can be justified for the respective examination.

CAUTION

Because of a technical error the acquisition is not stopped or there is an unintended release of radiation.

Risk of unnecessary radiation exposure!

- ◆ In the event of unwanted radiation, press the nearest EMERGENCY STOP button.

CAUTION

Failure of radiation indicator.

Danger of unintended radiation!

- ◆ If a radiation indicator fails contact Siemens Customer Service.

CAUTION

Additional objects in the beam path.

Risk of unnecessary radiation exposure!

- ◆ Note that additional objects in the beam path increase scatter radiation.

2 Safety

Please be aware that certain materials in the X-ray beam (e.g. parts of an operating table) may impair the X-ray image due to imaging of contours and inclusions in these materials. In certain rare cases, this may lead to incorrect diagnosis. This material may also result in higher radiation exposure.

Cios Fusion contains no patient positioning aids. Therefore, undertable radiation protection options are not included in the scope of supply. For example, overhanging lead rubber mats can be used to protect against scattered radiation underneath the patient table.

Deterministic radiation effect

As per IEC 60601-1-3, 5.2.4.5 (A.2), deterministic radiation injury is possible if a radiation dose delivered to an organ or tissue exceeds a value of 1 to 3 grays.

With typical applications and proper use of the C-arm, there is no reason to expect such radiation injuries. It is assumed that the maximum fluoroscopic time will not exceed 20 minutes, depending on the application, and that the point of skin penetration will be 30 cm (50 cm maximum) away from the detector input.

Example: An accumulated fluoroscopic time of 20 minutes and a skin penetration dose of 20 mGy/min yields a dose of 400 mGy.

Certain applications in the medical disciplines of vascular and cardiac surgery in particular require a significantly longer radiation exposure time of up to 60 minutes.

Example: The dose administered during individual treatments can exceed 1 Gy, while the average radiation dose in these procedures is normally below 500 mGy.

For procedures with a longer than expected radiation time, we recommend varying the beam direction (oblique beam through orbital and/or angular rotation) during the course of the procedure. In addition, the system reports the beam duration that has already been applied (regardless of applied dose).

The skin penetration dose for various operating modes and under standard operating conditions can be estimated using the "Dosimetric information" table (→ Page 251 *Dosimetric information*).

When changing the distance from the skin penetration point to the focus, please note that the skin dose decreases with the inverse square of the distance to the focus. This means that when the distance to the focus is halved, the skin penetration dose rate quadruples.

CAUTION

Remaining in hazardous area.



Danger of inadvertent radiation!

- ◆ During radiation release, make sure you are in a location subject to as little radiation as possible.

Location and size of the relevant operating areas

All examination types may be performed in the illustrated operating areas of the X-ray system.

X-ray tube assembly at the bottom

2 Safety

Maximum scatter radiation in operating area Scatter radiation in the main operating area according to EN 60601-1-3

Height above the floor [cm]	Measurement A1 [mGy/h]	Measurement A2 [mGy/h]	Measurement B1 [mGy/h]	Measurement B2 [mGy/h]
10	2.10	2.35	0.14	0.67
20	2.51	2.84	0.13	0.68
30	3.25	3.34	0.12	0.68
40	4.60	3.66	0.12	0.68
50	5.16	4.24	0.13	0.69
60	5.71	5.14	0.11	0.68
70	6.79	6.07	0.14	0.67
80	8.24	7.30	0.15	0.65
90	10.07	8.64	0.19	0.59
100	11.74	9.33	0.18	0.53
110	10.69	7.82	0.16	0.45
120	8.97	5.40	0.15	0.41
130	7.41	4.56	0.15	0.40
140	3.24	2.51	0.14	0.39
150	1.47	1.50	0.14	0.36
160	0.71	1.02	0.13	0.34
170	0.47	0.72	0.13	0.30
180	0.35	0.55	0.11	0.27
190	0.28	0.44	0.11	0.23
200	0.24	0.30	0.11	0.20

■ Measurement A1: Operating area A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position down, with scattered radiation grid

■ Measurement A2: Operating area A2


Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position down, with scattered radiation grid

■ Measurement B1: Operating area B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position down, with scattered radiation grid

■ Measurement B2: Operating area B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position down, with scattered radiation grid

(1) Measurement A

(2) Measurement B

Measurement A

■ Measurement A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement A2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Measurement B



■ Measurement B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

X-ray tube assembly at the top

Maximum scatter radiation in operating area

Scatter radiation in the main operating area according to EN 60601-1-3

Height above the floor [cm]	Measurement A1 [mGy/h]	Measurement A2 [mGy/h]	Measurement B1 [mGy/h]	Measurement B2 [mGy/h]
10	0.47	0.71	0.08	0.30
20	0.78	0.93	0.09	0.33
30	1.31	1.32	0.10	0.37
40	2.73	2.14	0.10	0.39
50	6.79	3.83	0.08	0.40
60	8.18	4.65	0.08	0.40
70	10.09	6.77	0.07	0.44
80	11.51	8.99	0.08	0.51
90	10.33	8.93	0.14	0.58
100	8.64	7.95	0.18	0.63
110	7.21	6.68	0.21	0.66
120	6.13	5.71	0.22	0.67
130	5.54	4.73	0.19	0.66
140	4.99	4.03	0.24	0.65
150	3.67	3.69	0.24	0.65
160	2.71	3.10	0.22	0.63
170	1.79	2.47	0.19	0.62
180	1.49	1.98	0.22	0.58
190	1.30	1.52	0.24	0.55
200	1.09	1.21	0.25	0.53

■ Measurement A1: Operating area A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position up, with scattered radiation grid

■ Measurement A2: Operating area A2


Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position up, with scattered radiation grid

■ Measurement B1: Operating area B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position up, with scattered radiation grid

■ Measurement B2: Operating area B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position up, with scattered radiation grid

(1) Measurement A

(2) Measurement B

Measurement A

■ Measurement A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement A2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Measurement B



■ Measurement B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Detector, lateral 40°

Maximum scatter radiation in
operating area

Scatter radiation in the main operating area according to EN 60601-1-3

2 Safety

Height above the floor [cm]	Measurement A1 [mGy/h]	Measurement A2 [mGy/h]	Measurement B1 [mGy/h]	Measurement B2 [mGy/h]
10	0.98	1.49	0.22	0.83
20	1.06	1.84	0.15	0.88
30	1.18	2.38	0.24	0.95
40	1.32	3.08	0.26	1.00
50	1.48	3.89	0.21	1.06
60	1.64	5.03	0.18	1.09
70	1.74	6.74	0.35	1.11
80	1.81	8.31	0.37	1.10
90	2.10	10.76	0.43	1.14
100	2.65	15.12	0.43	1.11
110	2.86	20.46	0.42	0.98
120	2.02	18.73	0.40	0.89
130	1.04	15.51	0.37	0.81
140	0.75	8.58	0.35	0.73
150	0.60	3.55	0.32	0.65
160	0.59	2.30	0.29	0.59
170	0.59	1.61	0.28	0.50
180	0.54	1.14	0.25	0.43
190	0.49	0.84	0.22	0.37
200	0.43	0.62	0.19	0.31

■ Measurement A1: Operating area A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical -40°, X-ray tube position lateral -40°, with scattered radiation grid

■ Measurement A2: Operating area A2


Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical -40°, X-ray tube position lateral -40°, with scattered radiation grid

■ Measurement B1: Operating area B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical -40°, X-ray tube position lateral -40°, with scattered radiation grid

■ Measurement B2: Operating area B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical -40°, X-ray tube position lateral -40°, with scattered radiation grid

(1) Measurement A

(2) Measurement B

Measurement A

■ Measurement A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement A2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Measurement B

■ Measurement B1



Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

2 Safety

X-ray tube assembly lateral 45°

Maximum scatter radiation in operating area

Scatter radiation in the main operating area according to EN 60601-1-3

Height above the floor [cm]	Measurement A1 [mGy/h]	Measurement A2 [mGy/h]	Measurement B1 [mGy/h]	Measurement B2 [mGy/h]
10	0.94	1.10	0.06	0.29
20	1.15	1.23	0.07	0.30
30	1.21	1.47	0.07	0.31
40	1.41	1.75	0.07	0.32
50	1.69	2.34	0.06	0.32
60	2.14	3.69	0.06	0.33
70	2.69	5.42	0.08	0.36
80	2.95	7.01	0.11	0.40
90	2.85	8.70	0.14	0.43
100	3.12	10.11	0.14	0.45
110	3.08	12.51	0.14	0.46
120	2.98	16.91	0.13	0.45
130	2.68	20.58	0.12	0.41
140	2.45	18.52	0.12	0.37
150	2.16	11.72	0.11	0.35
160	1.85	5.87	0.11	0.32
170	1.57	3.99	0.11	0.30
180	1.27	2.81	0.10	0.28
190	1.05	1.96	0.10	0.25
200	0.83	1.38	0.10	0.23

■ Measurement A1: Operating area A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position lateral 45°, with scattered radiation grid

■ Measurement A2: Operating area A2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position lateral 45°, with scattered radiation grid


■ Measurement B1: Operating area B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position lateral 45°, with scattered radiation grid

2 Safety

■ Measurement B2: Operating area B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm vertical, X-ray tube position lateral 45°, with scattered radiation grid

(1) Measurement A

(2) Measurement B

Measurement A

■ Measurement A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement A2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Measurement B



■ Measurement B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

X-ray tube assembly horizontal

Maximum scatter radiation in
operating area

Scatter radiation in the main operating area according to EN 60601-1-3

2 Safety

Height above the floor [cm]	Measurement A1 [mGy/h]	Measurement A2 [mGy/h]	Measurement B1 [mGy/h]	Measurement B2 [mGy/h]
10	1.27	1.62	0.09	0.35
20	1.37	2.00	0.10	0.36
30	1.50	2.65	0.09	0.37
40	1.64	3.53	0.06	0.40
50	1.86	4.85	0.05	0.43
60	2.16	6.72	0.06	0.46
70	2.81	9.11	0.05	0.50
80	3.16	11.25	0.05	0.52
90	2.89	11.88	0.07	0.53
100	2.51	10.45	0.07	0.52
110	2.30	8.03	0.07	0.49
120	2.10	5.86	0.09	0.45
130	1.72	4.14	0.11	0.41
140	1.45	3.02	0.09	0.37
150	1.27	2.22	0.11	0.33
160	1.10	1.68	0.10	0.29
170	0.95	1.28	0.10	0.28
180	0.83	0.99	0.10	0.28
190	0.71	0.78	0.10	0.25
200	0.62	0.59	0.09	0.24

■ Measurement A1: Operating area A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm horizontal, detector position lateral 90°, with scattered radiation grid

■ Measurement A2: Operating area A2


Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm horizontal, detector position lateral 90°, with scattered radiation grid

■ Measurement B1: Operating area B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm horizontal, detector position lateral 90°, with scattered radiation grid

■ Measurement B2: Operating area B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor; C-arm horizontal, detector position lateral 90°, with scattered radiation grid

(1) Measurement A

(2) Measurement B

Measurement A

■ Measurement A1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement A2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Measurement B

■ Measurement B1

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

■ Measurement B2

Fluoroscopy 110 kV, 3.0 mA continuous, 18 cm x 18 cm at the entrance of image receptor

Radiation interruption for all operating modes

The hand switches and footswitches are designed as pushbutton switches. Radiation is interrupted in fluoroscopy when releasing the corresponding operating element or, in other operating modes, after acquiring the stored image.

2.3 Equipment Safety

2.3.1 Mechanical safety

CAUTION

System component error.

An error on one system component can reduce the safety of the entire system.

- ◆ Terminate system operation.
- ◆ Contact Siemens Customer Service.

CAUTION

Image intensifier housing or detector cover damaged.

Toxic contamination!

- ◆ If the image intensifier housing or the detector is damaged, do not continue using it.
- ◆ Contact Siemens Customer Service.

CAUTION

Spindle nut broken or worn.

Risk of crushing!

- ◆ Vertical lift is blocked if the mechanical safeguard was triggered.
Contact Siemens Customer Service to replace the mechanical safeguard.

2.3.2 Positioning the C-arm

In case of improper handling of the Cios Fusion, the mobility of the C-arm may lead to collisions of the image intensifier/flat detector and the single tank with the patient and the patient table.

Brakes

Make sure the brakes are applied after adjusting the C-arm position.

Transport

When moving or transporting the C-arm system please take special care that the system parts do not collide with an obstacle. This could also result in accidental radiation release or an impairment of image quality under certain circumstances.

CAUTION

Transport damage.

Risk of injury!

- ◆ If repairs are not done by a professional, injuries to the user, patient or third parties and consequential product damage cannot be ruled out.

2.3.3 Installation, repair

Modifications of or additions to the product must be made in accordance with the legal regulations and generally accepted engineering standards.

As the manufacturer, Siemens cannot accept responsibility for the safety features and for the reliability and performance of the equipment if:

- The product is used in a manner other than that specified in the Operator Manual
- Installation, upgrades, readjustments, modifications or repairs are performed by personnel not contracted and authorized by Siemens
- Components affecting safe operation of the product are not replaced by original spare parts in the event of a malfunction
- The electrical wiring in the room containing the system does not meet the specifications of DIN VDE 0107 or the corresponding local regulations

If desired, we will provide the technical documentation for the product. However, this does not imply authorization to undertake repairs.

We cannot be held responsible for repairs made without our express written approval.

When any work is performed on the product, we recommend that you obtain a certificate indicating the nature and scope of the work performed. The certificate should include any changes in rated parameters or operating ranges as well as the date, the name of the company and a signature.

2.3.4 Original accessories

For safety reasons, only approved original accessories or non-Siemens accessories approved by Siemens Healthcare GmbH, may be used for this product.

The operator is liable for any risks associated with the use of accessories not approved by Siemens.

CAUTION

Inappropriate accessories.

The use of accessories that do not comply with the safety requirements of this equipment can reduce the safety of the entire system.

- ◆ Use only original Siemens accessories or accessories approved by Siemens.

2.3.5 Combination with other systems

To ensure the required safety, only products/components expressly approved by Siemens Healthcare GmbH may be used in combination with this system.

Please note that changes to the system may be carried out only with the express authorization of Siemens Healthcare GmbH.

Additional components placed into the beam path (e.g. positioning aids) will attenuate radiation and can degrade image quality.

2.3.6 Attachment of dedicated options

The attachment of certain (dedicated) options is permitted only if the following conditions are met:

General safety requirements

The use of accessories that do not comply with the relevant safety requirements of this system can result in a reduced safety level of the combined system.

When choosing accessories, the following aspects must be considered in particular:

- Use of accessories close to the patient.
- Proof that the accessories have been safety tested according to the applicable IEC 60601-1 guideline and/or the IEC 60601-1-1 harmonized national standard.

CAUTION

Changing the overall center of gravity.

Risk of crushing!

- ◆ Only products/components approved by Siemens may be installed.

Attenuation equivalent

According to IEC 60601-1-3, inadequate attenuation of the X-ray beam by materials between the patient and image receptor must be avoided.

Documented proof by the manufacturer is recommended.

Remove any auxiliary devices located in the beam path for calibration or adjustment of the dedicated options before operating the Cios Fusion.

Image quality can be impaired by placing materials directly in front of the image intensifier/detector, or the applied dose is increased by the automatic adjustment. Additional objects in the beam path may result in increased scattered radiation.

Weight counterbalance

CAUTION

Changing the weight distribution

Risk of crushing!

- ◆ Adding weight to the image intensifier/detector or single tank side means a loss in counterweight and can result in unintended movement of the C-arm.

Users must be alerted to the loss of counterbalance by a warning label. The responsibility for affixing the corresponding warning label lies with the company that attaches the dedicated option to the C-arm.

Image quality

The attachment of a dedicated option must not affect image quality (impairment of the follow-up).

After maintenance or service work, the correct function of the non-Siemens system on the detector must be tested.

Electrical safety

EN 60601-1, Section 3, "Protection against electric shock hazard" must be complied with.

Electromagnetic compatibility

EN 60601-1-2 must be observed in order to comply with the limit values for electromagnetic compatibility.

Attaching navigation systems

CAUTION

Additional weight after installing components for navigation.

Risk of injury!

- ◆ Note that the C-arm is no longer balanced after additional components are attached. It can move on its own after brakes are released.

CAUTION

Heavy weight of external navigation components.

Risk of injury!

- ◆ Note that stability in accordance with DIN EN/UL standards is guaranteed only if the external system on the image intensifier/detector (2D navigation) does not exceed the calculated maximum weight.

⚠ CAUTION

Overloading the C-arm.

Risk of injury!

- ◆ Note that the braking action is guaranteed only if the external system on the image intensifier/detector (2D navigation) does not exceed the calculated maximum weight.

⚠ CAUTION

Sharp edges.

Risk of injury!

- ◆ Be careful of any sharp edges on the externally attached navigation system.

⚠ CAUTION

Incorrectly attached navigation system.

Risk of injury!

- ◆ Before navigation-assisted operation, the user must ensure sufficient accuracy of the navigation system.

⚠ CAUTION

Insufficient navigation accuracy due to contaminated or defective detector interface.

Risk of injury!

- ◆ Prior to using the navigation system check the required cleanliness and functionality.

⚠ CAUTION

Use of insufficiently electrically-protected external navigation components.

Risk of injury!

- ◆ Use only accessories corresponding to the CE Declaration of Conformity.

Additional safety information

- To avoid thermal overloading of components and short circuits, EN 60601-1 Section 7 and, if appropriate, UL 60601-1 must be complied with.
- Connecting external loads to the power supply of the C-arm system is not permitted.
- We recommend that users in the EU have the relevant manufacturer of the accessory operated by you confirm the CE Declaration of Conformity according to Appendix II, MDD (Directive 93/42/EEC dated June 14, 1993), Directive 2011/65/EU dated June 8, 2011 on the restricted use of certain hazardous substances in electrical and electronic devices, and the Declaration of Compatibility according to Article 12, MDD. In countries outside the EU the relevant national regulations must be observed.

The product liability and warranty are restricted or expire if the above listed conditions and limit values are not complied with when attaching accessories.

For non-Siemens options we generally accept no liability.

2.3.7 Disposal

- If you want to remove the product from service, take into consideration that public legal directives may contain special regulations regarding disposal of this equipment. In order to ensure that these legal regulations are complied with and to avoid potential environmental hazards which may be caused by the disposal of your system, please consult Siemens Customer Service.
- Batteries and packaging material must be disposed of in an environmentally safe manner according to national regulations.
- The sterile single-use covers must be disposed of in accordance with national regulations or the rules of the hospital.
- For further information about the disposal of the product, please refer to our service documents.

2 Safety

3 System Description

3.1	Device Description	51
3.1.1	System overview	51
	Options	51
3.1.2	C-arm system	53
	Flat detector with grid	53
	Single tank	54
3.1.3	Monitor trolley	55
	UPS display field	56
3.1.4	Control elements	56
	Control unit on the C-arm system	57
	Control panel on the monitor trolley	58
	Hand switch	59
	Footswitch	59
3.1.5	Operating modes	60
	Fluoroscopy	60
	Single image	61
	Subtraction/Roadmap (optional)	61
3.1.6	Description of image processing	62
3.1.7	Image display	62
3.2	System Operation	63
3.2.1	Start-up	63
	Connecting the C-arm system with the monitor trolley	63
	Connecting the footswitch	64
	Establishing the equipotential bonding connection	65
	Establish the power line connection	65
	Switching on the Cios Fusion	66
	Functions during system start-up	66
3.2.2	Move C-arm	67
	Operating the brakes	69
	Lifting and lowering the C-arm	69
	Moving the C-arm horizontally	71
	Swivelling the C-arm	72
	Angulating the C-arm	72
	Orbital movement of the C-arm	73
3.2.3	Control panel	74
	Exposure parameters	74
	Task cards	75
	Control elements	76
	Submenus	77
	Keyboard mode	78
3.2.4	Operating the hand switch	79
	Radiation release	79
	Image storage	79
	Storage	80
3.2.5	Operating the footswitch	80
	Radiation release	80
	Image storage	80


3 System Description

3.2.6	Switching off	81
	Switching off the system completely	81
	Switching off the C-arm system only	82
	Disconnecting Cios Fusion from line power	82
	Reactivating the system	83
3.2.7	Transporting the C-arm system	83
	Preparing the C-arm system	83
	Drive operation and steering of the C-arm system	84
	Braking the C-arm system	85
	C-arm system park position	86
3.2.8	Transporting the monitor trolley	86
	Preparing the monitor trolley	86
	Moving the monitor trolley	87
	Locking the monitor trolley	87
	Setting the direction of the monitor trolley	88
	Monitor trolley park position	88
	Park position of the monitors with height-adjustable monitor column (optional)	88
3.3	Emergency Situations	90
3.3.1	Function of the Emergency Stop button	90
3.3.2	Behavior in case of power failure	91
	Restarting	91
3.3.3	Restarting the Cios Fusion	91
3.3.4	Restart via reset	92
	Resetting	92
3.3.5	Procedure for cardiopulmonary resuscitation (CPR)	93

3.1 Device Description

3.1.1 System overview

The Cios Fusion consists of a C-arm system and a monitor trolley.

- (1) C-arm system with flat detector (20 cm x 20 cm or 30 cm x 30 cm depending on the configuration) and stationary anode tube with generator
- (2) Monitor trolley with touch screen, mouse, USB port, two rotatable TFT displays, DVD drive, and memory for 150,000 images

Options

The following options are available for the Cios Fusion:

Monitor trolley

- Height-adjustable and foldable (Park position) monitors
- 2 TFT Premium high-brightness color displays (alternative to standard TFT high-brightness color displays)
- Wireless network connection (WLAN)
- Sony UP-D898MD thermal printer
Economical, light-weight, compact printer for A6 b/w printing on thermal paper
- SONY UP-971AD, UP-991AD thermal printer
Economical, light-weight printer for A4 b/w printing on thermal paper (UP-971AD, UP-991AD) and transparent film (UP-991AD)
- Audio package


C-arm

- 30x30 detector, 20x20 detector
- Position memory

3 System Description

- Integrated detector laser light localizer
- Single-tank laser light localizer
- Dose measurement chamber for dose area product/air kerma (valid for IEC 60601-1:2005 only)
- Control elements**
 - Remote control unit
 - Additional holder for remote control unit
 - Multifunctional footswitch (with extended functionality)
- Applications**
 - Vascular Software
 - Vascular Software Premium
 - Stenosis quantification
 - Live graphic overlay (drawing lines and displaying on live image)
 - Radiation lock code
 - 2D measuring function (for measuring angles and distances)
 - DICOM Worklist
 - DICOM Send
 - DICOM Print
 - DICOM Query/Retrieve
 - HIPAA (Health Insurance Portability and Accountability Act)
 - LithoReport
 - CD/DVD Recording (storage of fluoro sequences on CD/DVD in MPEG4 format)
- Interfaces**
 - DVI video splitter (Monitors A and B)
 - SmartView - HD video manager
 - CAT DVI
 - Interface for injection pump triggering
 - NaviLink 2D (integrated digital navigation interface for lossless transfer of 2D image data to a navigation system)
 - Litho option (for connecting a MODULARIS Variostar ESWL system)
- Accessories**
 - Grounding cable
 - DHHS spacer
 - Sterile covers for the detector and single tank
 - Sterile cover with plastic brackets for the C-arm
 - Sterile cover for the remote control unit
 - Non-sterile cover for (standard/multifunctional) footswitch

3.1.2 C-arm system


- (1) C-arm
- (2) Handles for manual movement/positioning of the C-arm
- (3) Control unit
- (4) Horizontal support arm
- (5) Handles on both sides with a holder for the hand switch, as well as a steering handle (center) for moving/transporting and braking the C-arm system
- (6) Lifting column
- (7) Cios Fusion Electronics unit
- (8) Footswitch (standard)
- (9) Flat detector with grid
- (10) Single tank with X-ray tube unit and integrated collimator
- (11) Wheels with cable deflectors

Flat detector with grid

Compared with the classic image intensifier, the flat detector has a much higher dynamic range. This enables an expansion of the application range, for example, to soft tissue imaging.


The grid attached to the flat detector further increases image quality.

3 System Description


- (1) Flat detector
- (2) Handle for manual movement/positioning of the C-arm
- (3) Removable grid

Single tank

- (1) Exit window for optional laser light localizer
- (2) Integrated collimator system

3.1.3 Monitor trolley

- (1) Radiation display
- (2) Monitors can be rotated by 180°
(optional: folding monitors)
- (3) Control panel
- (4) Handles: for neat cable storage
- (5) Drawer for storage of CDs/DVDs, Quick Guide, etc.
- (6) DVD R/W drive
- (7) Wheels with cable deflectors
- (8) Monitor column
(optional: height adjustable)
- (9) Storage for mouse, MP3 player;
Cios Fusion on/off switch, connections for USB, mouse
- (10) Connection for WLAN, LAN, Reset key
- (11) Shelf for printer
- (12) UPS window (Charge state)
- (13) Central locking brake


3 System Description

UPS display field

(1) UPS window

The UPS display field on the front of the monitor trolley shows the charge state of the batteries and the operating status.

(1) LED "Green": Input voltage (line power operation)

(2) LED "Yellow": Battery operation

(3) LED "Red": Alarm, malfunction

(4) Bar display of utilization ratio


(5) Bar display of battery charge status

When the red LED is lit, the Cios Fusion is no longer protected against power outages by the UPS. The power plug may not be pulled until the imaging system is shut down. If this is the case, notify Siemens Service.

3.1.4 Control elements

The following graphic is an overview of the Cios Fusion control elements. Some of the control elements presented here are available as options.

- (1) Control unit on the C-arm
- (2) Hand switch
- (3) Control panel on the monitor trolley
- (4) Optional: Remote control unit with cable
- (5) Standard footswitch
(optional: multifunctional footswitch)


Use only your fingers to operate the touch screens. Do not use hard or sharp objects.

Control unit on the C-arm system

The control unit for performing your examinations is located on the C-arm system.

3 System Description

- (1) Power supply indicator (LED)
- (2) Radiation indicator
- (3) EMERGENCY STOP
- (4) Vertical movement of the C-arm
- (5) Control panel


The control panel on the C-arm control unit offers the same functions and shows the same user interface as the control panel on the monitor trolley.

The control panel is not an imaging device as prescribed in DIN 6868-157. Any images displayed are not suitable for diagnostic purposes.

Control panel on the monitor trolley

The monitor trolley control panel is used to operate applications for preparing (e.g. entering patient data) and evaluating examinations. The control panel is integrated flush into the cover plate of the monitor trolley.

(1) Control panel on the monitor trolley

The control panel is not an imaging device as prescribed in DIN 6868-157. Any images displayed are not suitable for diagnostic purposes.

There is also a mouse, which is used as usual to execute functions on the left monitor.

Hand switch

- (1) Yellow button for radiation release
- (2) Holder for hand switch
- (3) Key for saving images

Footswitch

The footswitch pedals assignment can be configured.

3 System Description

Configuration for foot pedal assignment

- (1) Radiation release in preselected operating mode (standard setting)

Exception: Single image when fluoroscopy operating mode is preselected

- (2) Radiation release for fluoroscopy (Standard setting)

In addition, the image saving function can be allocated to a pedal if desired. (Standard footswitch only).

3.1.5 Operating modes

The following presents the operating modes of the Cios Fusion, their properties, and assigned system settings.

Fluoroscopy

For **Fluoroscopy** operating mode, you can choose between several examination sets with different characteristic curves for fluoroscopy to determine the radiation parameters for fluoroscopy. Every application allows you to choose between different optimized programs.

The pulse duration is generally between 7 and 40 milliseconds. According to the level of noise reduction, many different fluoroscopic images can be integrated.

The standard setting of the Cios Fusion after startup is established in the configured standard application, and in general is **Fluoroscopy** (→ Page 219 *Examination settings (PEX Editor)*).

System settings

Exposure factors and system control units, including the way in which the automatic setting is controlled:

- 1 k² matrix (1.5 k² matrix for 30 cm x 30 cm detector);
- Frame rate usually 0.5 to 30 frames/s;
- Image integration (as a function of the k factor set), i.e., a number of k exposures are integrated into one image by sliding averaging; the k factor can be selected between k = 1 ("OFF" setting) and k = 16 (32 for LIH) and can be assigned to an examination set and stored.

Typical clinical procedure

Fracture reposition of the distal upper extremity (e.g. distal forearm fracture) in the plaster room of an emergency outpatient clinic where, under fluoroscopy, the fracture elements are reduced by extension, fixed temporarily in the best possible position and then fixed permanently by applying a plaster cast.

Single image

The **Single image** operating mode provides an electronic instant image of the patient on the monitor. It is recommended for final exposures. The exposure time depends on the pulse width set.

System settings

Exposure factors and system control units, including the way in which the automatic setting is controlled:

- 1 k² matrix (1.5 k² matrix for 30 cm x 30 cm detector);
- X-ray pulse with a width of 7 ms up to a maximum of approx. 1800 ms.

Typical clinical procedure

Final follow-up exposure of a fracture reposition of the distal upper extremity (see above).

Subtraction/Roadmap (optional)

The subtraction memory option allows you to perform a digital subtraction angiography and simultaneously display the unsubtracted angiogram on the second monitor. Subtraction technique allows hemodynamic display as well as display of the maximum vascular filling and Roadmap. The Roadmapping features can also be used for other procedures.

System settings

Exposure factors and system control units, including the way in which the automatic setting is controlled:

- 1 k² matrix (1.5 k² matrix for 30 cm x 30 cm detector);
- Acquisition rate is freely selectable, usually 4 to 15 frames/s;
- Image integration (as a function of the k factor set), i.e. a number of k exposures are integrated into one image; the k factor can be set between k = 1 and k = 16 (32 for LIH) by an authorized technician.

Typical clinical procedure

- Display of an arterial vessel for localizing vascular stenoses with injection of a contrast agent to enable the contrast-enhanced display of the vascular filling (subtraction of the native image (mask) from the contrast-enhanced image).
- Alternative to native display, subsequent inversion of the displayed image allows you to display a catheter introduced into the vessel path using the Roadmap function.

For information on performing an examination with Subtraction or Roadmap, refer to (+ Page 133 *Performing special examinations*).

3.1.6 Description of image processing

All image processing of the modality is executed by the "Imaging chain" component, which is part of the image system software. The imaging chain is divided into hardware and software modules. It is responsible for correct parameterization of the hardware modules (e.g. the flat detector or the image intensifier group). From there, the raw image data from the image converter is transferred via exclusive image transfer interfaces to the PC hardware using streaming technology. On the PC, the image data is then analyzed and processed in real time, and an initial adjustment of the brightness and contrast is performed. Finally, the image data is transferred to the image system software where it is stored. At the same time, the image data is processed for optimum medical evaluation for the application in question (including local filtering, subtraction, brightness and contrast processing, geometric operations, "Peak opacification") and displayed without delay on one or more monitors. The export for the preview, on touch displays for example, is done in the same way (if supported by the system).

If a DICOM report is needed, the image data can be reimported into the process. All required processing steps are then performed again and the image data re-exported to the imaging system software. From there they are sent to a (DICOM) archive along with the patient data.

Note that the exported images are not transmitted as raw data, but include all of the applied processing steps. The original condition cannot be restored (please do not use these images for primary diagnosis!).


3.1.7 Image display

The imaging system uses two monitors for image display.


- The left monitor is used to display live images as well as preoperative and postoperative images.
- The right monitor is used to display the reference images.

Reference images may be:

- preoperative images, incl. images of other modalities (e.g. CT, MR)
- images of current examinations
- images from previous examinations

3.2.1 Start-up


The C-arm system is connected to the monitor trolley with a cable.

Tripping over cable.

Risk of falling.

- 63


3 System Description

- 1 Plug the central plug into the socket on the left side of the C-arm system, when viewed from the C-arm. The switch must always be in the vertical position.

The monitor trolley may only be connected to the corresponding C-arm system. If the monitor trolley is connected to the wrong C-arm system, an error message is displayed indicating that the C-arm and monitor trolley are not a suitable pairing. Upon acknowledging the message, the text "Incorrect C-arm" indicating the status is displayed at the lower left of the imaging system monitor.


- 2 Turn the central plug switch to the right until it audibly clicks into place.

The monitor trolley is connected to the C-arm system.

Connecting the footswitch

For the release of radiation with the footswitch, it must be connected to the C-arm system.


- ◆ Plug the footswitch cable into the socket labeled with the footswitch symbol on the front of the C-arm system.

Establishing the equipotential bonding connection

The Cios Fusion can be connected to a protective ground terminal via the equipotential bonding connector on the C-arm system. This will ensure that the Cios Fusion has the same electrical potential as other units connected to the same protective ground terminal.

When performing cardiac examinations or examinations of the open skull, an additional grounding cable according to DIN 57107/VDE107 must be routed in rooms of Application Group 2.

- ◆ Clamp the equipotential bonding cable to the socket indicated by the symbol on the C-arm system front connection panel (arrow) and to an equipotential bonding point in the vicinity of the patient.

Equipotential bonding is established.

Establish the power line connection

The Cios Fusion operates using a line power connection cable from the monitor trolley to a properly grounded socket.

3 System Description

CAUTION

Cable damage.

Risk of electrical shock!

- ◆ Check the power cable. Do not use the device if the power cable is damaged.

- 1 Plug the monitor trolley power plug into the appropriate socket.

The line power connection is established.

- 2 Position the unit so that it is easy to disconnect from the line power.

Switching on the Cios Fusion

- ◆ Press the right **ON** button on the monitor trolley.

The button that can be used to make a change is always illuminated. Thus, when the system is shut off, the ON button is illuminated.

Functions during system start-up

The Cios Fusion is switched on.

The system automatically runs a self-test.

On the control panel, all LEDs light for approx. 2 seconds as a function test and then switch off automatically.

Check whether all LEDs light correctly. If not, notify Siemens Customer Service as there could be a malfunction.

An audible signal sounds while the system starts up.

If an audible signal does not sound, notify Siemens Customer Service; there is a possibility that audible signals will not play back correctly.

An automatic self-test is performed for the C-arm control unit buttons.

Do not press any keys (such as the lift key) during booting; otherwise, an error message will be displayed indicating that a key was pressed during booting, and therefore the self-test of the keys did not function.

If error messages are displayed during the self-test, a new self-test is performed automatically. The self-test can be performed up to four times.

A failed self-test is indicated on the affected control panel by continuous flashing of all LED fields. Use of the buttons is not possible in this state; the affected control panel is then shut down automatically. In this case, Siemens Customer Service must be notified.

Radiation stand-by

On the left monitor, radiation stand-by is indicated by a symbol in the status bar (lower left). On the control panel, radiation stand-by is indicated by a green horizontal bar in the exposure parameters area (top left).

This icon appears on the monitor if radiation release is blocked. On the control panel, the horizontal bar is grey. You need to enter the password to unlock the radiation release (→ Page 101 *Unlocking/locking radiation release*).

Radiation cannot be released if the monitor trolley is operated separately. On the control panel, the horizontal bar is white.

Essentially, the system is not ready for radiation unless a patient is registered.

If the C-arm system is connected to a monitor trolley that is already switched on, it normally becomes available after 120 s. It is then in the same operating state as it was prior to being disconnected from the monitor trolley.

Flat detector calibration message

If a message is displayed on the monitor indicating that flat detector calibration is due:

- 1 Perform detector calibration.
- 2 Confirm it to continue start up.

All images acquired with an uncalibrated system are tagged "without valid calibration".

Presettings

The following functions are set when the Cios Fusion is started:

- Operating mode: the preset operating mode in the application configured as **Standard** (→ Page 219 *Examination settings (PEX Editor)*)
- Image mirroring: switched off if it is not switched on in the application configured as **Standard** (→ Page 219 *Examination settings (PEX Editor)*)
- Number of images (fluoro): the default value in the application configured as **Standard** (→ Page 219 *Examination settings (PEX Editor)*)
- Automatic dose rate control (ADR): Activated
- Zoom level: 0, if it is not switched on in the application configured as **Standard** (→ Page 225 *Editing applications*)
- Collimators: Full format

3.2.2 Move C-arm

The horizontal movement, swivel movement, angulation and orbital movement of the C-arm are performed manually. Use the handles on the C-arm or the flat detector.

Additionally, the C-arm can be adjusted in height by motor control.

3 System Description

WARNING

Manual or motorized movements.

Risk of crushing and collision for persons and objects in the vicinity of moving parts.

- ◆ During manual and motorized movements of the device, make sure you take into account any persons or objects in the way of movement, and control the movement to avoid collisions.
- ◆ Make sure that everyone is outside the hazardous zone.

CAUTION

Uncontrolled C-arm movements.

Risk of injury and crushing!

- ◆ Make sure the C-arm cable does not get caught in the horizontal support arm.

CAUTION

Activation of system movements without visual contact.

Risk of crushing!

- ◆ Use the system only if the patient and other persons can be observed during movements.

CAUTION

C-arm rotation.

Risk of crushing!

- ◆ Pay particular attention to crushing risks between moving system parts and the corresponding guide openings.

CAUTION

Movement due to unintended actuation of operating elements.

Risk of collision, risk of injury to the patient or operator, risk of damage to the device parts.

- ◆ Press the closest EMERGENCY STOP button if device movements do not stop.

Please note that the scale on the C-arm is intended for orientation only, not for measurement.

Operating the brakes

The Cios Fusion is equipped with mechanical brakes. The levers for releasing and locking the brakes for different directions of movement are marked with different colors. A graduated scale in the same colors for the corresponding directions of movement is located on the housing of the C-arm.


CAUTION

Brake failure.

Risk of crushing!

- ◆ Before beginning the examination, perform the daily function and safety checks.

Overview of C-arm brakes

- A: Horizontal brake
- B: Angular brake
- C: Swivel brake
- D: Orbital brake

Before moving the C-arm, the brake for the relevant direction of movement must be released.

Engage the corresponding brake lever after moving the C-arm.


Lifting and lowering the C-arm

You can lift and lower the C-arm by motor control using the arrow buttons on the C-arm system control unit.

3 System Description

The lifting column can be lowered from position 1 (highest level) to intermediate position 2, down to position 3 (lowest level).

- (1) Position 1 = highest level
- (2) Position 2 = intermediate position
- (3) Position 3 = lowest level

Lifting the C-arm

- ◆ Press the Up button on the C-arm system control unit.
The lifting column moves upwards.

Lowering the C-arm

The lifting column can be lowered to position 3.

- 1 Press the Down button on the C-arm system control unit.
The lifting column then moves to the intermediate position (2) and automatically stops there.
- 2 Press the Down button on the C-arm system control unit again.
The lifting column is lowered further.
An audible signal sounds.

Make sure there is nothing in the hazard zone (under the C-arm or travel frame) when lowering the C-arm; crushing hazard!

Lifting movement failure



If the lifting column cannot be moved the EMERGENCY STOP button is pressed and must be unlocked.


- 1 Turn the knob clockwise until it pops back out.
- 2 If the lifting column can still no longer be moved in any direction, then please contact Siemens Customer Service.

Moving the C-arm horizontally

You can move the support arm horizontally by up to 20 cm.



- 1 Open the brake lever marked in green for **horizontal movement** by turning it to position (1).
The brake is released.
- 2 Move the support arm to the desired position while observing the green scale.


- 3 Lock the brake lever by turning it back to position (2).
The brake is locked.

Swivelling the C-arm

You can move the C-arm in the horizontal plane $\pm 12^\circ$ about the lifting column.



- 1 Open the brake lever marked in orange for **horizontal swivel** by turning it to position (1).
The brake is released.
- 2 Swivel the C-arm to the required position.

- 3 Lock the brake lever by turning it back to position (2).
The brake is locked.


Angulating the C-arm

You can rotate the C-arm in the vertical plane $\pm 225^\circ$ about the horizontal support arm.

1 Open the brake lever marked in yellow for **angulation** by turning it to position (1).

The brake is released.



2 Rotate the C-arm to the required angulated position while observing the yellow scale on the support arm joint.

3 Lock the brake lever by turning it back to position (2).


The brake is locked.

Orbital movement of the C-arm

Starting from the basic position (0°), you can swivel the C-arm by up to +110° or -40° (150° in total).


3 System Description

- 1 Open the brake lever marked in blue for **orbital movement** by turning it to position (1).

The brake is released.

- 2 Swivel the C-arm to the required orbital position while observing the blue scale.


- 3 Lock the brake lever by turning it back to position (2).

The brake is locked.

3.2.3 Control panel

The control panels on the C-arm control unit, monitor trolley, and (optional) remote control unit are identical in terms of function and display.


Press the corresponding symbol button to operate the function.

- (1) Displays and operating elements for the exposure parameters
- (2) Task cards with operating elements for the current work step
- (3) **Close pat.** (close patient) button
- (4) Button to manually select keyboard mode

Exposure parameters

The left section of the control panel contains the operating and display elements for image acquisition. This area remains the same regardless of the task card selected.

(1) **Radiation stand-by indicator**

Green: ready, **white:** not ready, e.g., C-arm system not connected or patient not registered, **grayed:** Radiation release blocked (password entry required)

(2) **Display elements and position memory**

Current C-arm angles (for optional motorized operation and position memory), cumulative dose and radiation time, kV, mAs/mA, indication of grid status

(3) **Exposure preparation**

Collimator setting, laser light localizer, reset warning signal (radiation time)

(4) **Operating mode selection**

Fluoroscopy (Fluoro), Single image, Subtraction (Sub), Roadmap (Road)

(5) **Setting of radiation parameters**


Dose, Metal correction, Automatic dose control (ADR), zoom, noise filtering and motion detection, pulse frequency

Task cards

The task cards provide you access to all required operating elements in the current work step.

In the **PROCESSING** work step you are provided with operating elements for image postprocessing.

3 System Description

- (1) Basic operating elements for image display
- (2) Operating elements for special applications (subtraction, roadmap)
- (3) Operating elements for image postprocessing and -evaluation
- (4) Area for displaying submenus (→ Page 77 Submenus)
- (5) Image and information area

Changing task cards

- ◆ Press the tab of the required task card.

Control elements

You perform individual function and set the required parameters using the task card operating elements.

When you press a button and a dashed line appears, this means the requested function is being processed by the system. No additional operating steps are possible during this time for safety reasons.

Function value display: For some functions, the current set value is displayed in the operating element (example: Pulse rate 15 p/s).

Function inactive: Many functions can only be executed when specific conditions are met (for example, **Hold reference** only when the image display is on the right monitor. Functions that currently cannot be executed appear as inactive.

Performing the function

- ◆ Press the button for the required function.

For some functions the button can be pressed and held to continue performing the function (such as moving the collimator).

Activating/deactivating functions

- 1 Press the button in deactivated state to switch on the function.

- 2 Press the button in activated state to switch off the function.

Submenus

For functions requiring further input or additional settings, corresponding operating elements are available in a submenu.


Opening the submenu

- 1 Press the button for the required function.

The function is activated.

A submenu with function-specific operating elements is shown in the upper section of the task card.

- 2 Change the settings using the operating elements shown.

Close the submenu

- ◆ Press the Close button.

3 System Description

– or –

Press the button for the required function again.

Select a different submenu.

Keyboard mode

Keyboard mode is used to enter text and numbers.

- (1) Text field
- (2) Keyboard
- (3) Symbol button bar for opening functions

Opening keyboard mode

The system switches automatically to keyboard mode when functions requiring text input are opened (such as for patient registration). Keyboard mode can be opened manually when necessary (e.g., for service purposes).

- ◆ Press this button on the control panel (lower right).

A virtual keyboard appears on the control panel instead of the task cards.

Entering text

- ◆ Press the corresponding character buttons like a computer keyboard.

The entered characters are shown in the input screen on the monitor and in the white text field of the virtual keyboard.

For upper case letters and to enter special characters, press the **Shift** button and then the corresponding character.

To enter special characters briefly press in sequence: **Ctrl + Alt** or **Shift + Ctrl + Alt**

Ending keyboard mode

Keyboard mode ends automatically when the corresponding input screen on the monitor is closed. Depending on the function that is open, this can be done with the **Enter** button (such as comment text) or by opening the next work step (such as examination after patient registration).


- ◆ To terminate keyboard mode manually, press the Close button.
The task card previously shown is displayed again.

3.2.4 Operating the hand switch

The hand switch is used to remotely control radiation release and image storage.

The hand switch can be inserted into the holder provided on the C-arm system.

- (1) Radiation release key
- (2) Image storage key

Radiation release

The hand switch is used to release radiation in the preselected operating mode.

- ◆ Press the yellow radiation release key on the hand switch and hold it down while radiation is released.

The current radiation parameters are shown on the control units.

The currently generated image is displayed on the left monitor unless the **Hold reference** function is active.

Image storage

During radiation

- ◆ Press this key on the hand switch.
The image currently generated and displayed is saved.

After radiation

- ◆ Press this key on the hand switch.
Holding the key for < 2 seconds: Saves the image last recorded (LIH).
Holding the key for > 2 seconds: Saves the scene last recorded (LSH).

The images are stored in the local database. They are transferred from the left monitor to the right monitor unless the **Hold reference** function is active.

3 System Description

Storage

When not in use, the hand switch can be kept in one of the holders on either side of the C-arm system.

3.2.5 Operating the footswitch

The footswitch is used, for example, if both hands need to be free during the exposure.

The footswitch is also suitable for applications where fluids may land on the floor.

For operation of the multifunctional footswitch (optional) see (→ Page 291 *Multifunctional footswitch (optional)*).

Radiation release

In the default setting, the pedals of the footswitch are allocated as follows:

- The right pedal is always used to activate **Fluoroscopy** (Fluoro).
- The left pedal is used to activate the preselected operating mode.


Exception: If the **Fluoroscopy** (Fluoro) operating mode is preselected, the left pedal is allocated the **Single image** operating mode.

The left/right functionality of the two pedals can be reversed by Siemens Customer Service upon request.

Pedal allocation:

The current pedal allocation is shown on the monitor trolley left screen.

If the footswitch is defective or not connected, the symbol for the hand switch is displayed instead.

Releasing radiation

- ◆ Keep the left or right foot pedal pressed during radiation release.

The current radiation parameters are shown on the control units.

Image storage

One of the pedals can be configured with the "Save" function by Siemens Customer Service. In this case, in the default setting the pedals of the footswitch are allocated as follows:

- The right pedal is used to activate the preselected operating mode.
- The left pedal is used to activate image storage.

Saving during radiation

- ◆ Press the "Save" pedal on the footswitch.

The image currently generated and displayed is saved in the local database.

Saving after radiation

- ◆ Press the "Save" pedal on the footswitch.

Holding the pedal for < 2 seconds: Saves the image last recorded (LIH).

Holding the pedal for > 2 seconds: Saves the scene last recorded (LSH).

The images are stored in the local database. They are transferred from the left monitor to the right monitor unless the **Hold reference** function is active.

3.2.6 Switching off

Before disconnecting the Cios Fusion from line power, you must shut it down. During the shutdown procedure the imaging system is shut down before the system is switched off.

The shutdown procedure for the imaging system is complete when the stand-by indicator in the left monitor goes out.

CAUTION

Loss of image data.

Risk of unnecessary radiation exposure!

- ◆ Do not shut down the system while radiation is being released.

If you still need the monitor trolley (for postprocessing), you can disconnect the C-arm system while switched on. In this case, only the C-arm system is shut down.

Switching off the system completely

Ending active processes

- 1 Close the current patient with the corresponding button on the control panel.
- 2 Make sure data were not written to CD/DVD.
If necessary, wait for CD/DVD write processes to be completed.
- 3 If there is a disk in the CD/DVD, take it out of the drive.

Open network export jobs are retained during shutdown and are continued when the system is switched back on. Storage jobs for removable devices as well as print jobs are stopped and deleted.

3 System Description

Switching off

- ◆ Press the OFF button on the monitor trolley.
The C-arm system is immediately switched off.
The imaging system switches off after the computer is shut down.

The button that can be used to make a change is always illuminated. Thus, when the system is switched on, the OFF button is illuminated.

Switching off the C-arm system only

- ◆ Turn the central plug lever on the C-arm system to the left to unlock it. Unplug the connector.
The C-arm system is immediately switched off.
As long as a C-arm system is not connected, the acquisition functions on the imaging system are not available.

Disconnecting Cios Fusion from line power

Please note that after finishing an examination, the Cios Fusion must be shut down properly before it is disconnected from the power supply.

- 1 Switch off the Cios Fusion and wait for it to shut down.
- 2 Disconnect the monitor trolley power plug from the power outlet.

When the power plug is pulled, the UPS switches to battery mode (yellow LED illuminates).

Pull on the plug, not on the cable!

Releasing the connection cable


- 1 Turn the central plug lever on the C-arm system to the left to unlock it. Unplug the connector.
- 2 If attached, disconnect the equipotential bonding connection from the C-arm system.

Removing the remote control unit

- ◆ See (→ Page 294 *Removing the remote control unit and holder*).

Rolling up the footswitch cable

The footswitch is placed in the holder provided for storage and during transport.


- ◆ Roll up the footswitch cable (1) onto the cleat provided and place the footswitch into its holder (2).

When attaching the footswitch, please be careful not to kink the cable.

Reactivating the system

If the imaging system cannot be shut down properly after pressing the **OFF** button or if it no longer reacts to input, you must reactivate the Cios Fusion system as follows:

- 1 Press the **RESET** button on the right side of the monitor trolley.
All running processes are terminated and the Cios Fusion is shut down.
- 2 Switch the Cios Fusion on again and let it boot up completely.
Now you can either continue using the Cios Fusion or shut it down.

If the Cios Fusion is not fully operational despite the reset, please notify Siemens Customer Service.

3.2.7 Transporting the C-arm system


The C-arm system is equipped with 4 wheels for easy steering in any direction. The C-arm system can be locked in place with the steering/brake lever.

Preparing the C-arm system

Prior to transport, the C-arm system must be set to the transport position.

- 1 Release all brakes of the C-arm.
- 2 Set the C-arm to the transport position shown in the drawing.
Angulation 0°/Orbital position 0°.
- 3 Move the lifting column to its lowest position.
- 4 Move the horizontal support arm all the way back.

3 System Description

Drive operation and steering of the C-arm system

With the steering lever, you operate the breaks of the C-arm wheels and steer the C-arm in the required direction.

The steering lever can be locked into 3 different positions. One is for movement straight ahead, the others for transverse travel to the right or left.


CAUTION

Weight of the device.

Risk of crushing!

- ♦ Fast braking can be impacted by the weight of the device. Pay attention to your speed as well as the floor covering, slope, and uneven areas.


- 1 Press the release button and pull the lever up (indicator = green).


Drive operation for the C-arm system is now possible.

- 2 Turn the steering lever in the required direction.
- 3 Move the C-arm system by hand.

You have the following options for steering the C-arm system:

**■ Motion straight ahead / pivoting in place**

The brake is released, the steering lever is pointing straight up.
The C-arm system can now be driven straight ahead or pivoted in place.

**■ Diagonal motion**

The brake is released, the steering lever is angled to the left or to the right.
The C-arm system can now be driven diagonally.

**■ Straight left or straight right motion**

The brake is released, the steering lever is pointing straight to the left or the right.
The C-arm system can now be driven straight to the left or straight to the right.

Avoid transporting the C-arm system on inclined surfaces $> 5^\circ$ or it may slip or tip over.


When transporting the C-arm system make sure there are no obstructions on the floor.

Braking the C-arm system

To securely brake the C-arm system, particularly on tilted or uneven planes and surfaces, proceed as follows:

- ◆ Keep the unlocking button pressed and lift the steering lever.

The red indicator appears.

The C-arm system brakes are now locked. However, even with the brake applied, C-arm movements are still possible.

3 System Description

C-arm system park position

CAUTION


Insufficient braking force.

Risk of crushing!

- ◆ When parking the C-arm system on an inclined surface, lock the brakes on the C-arm system and turn the wheels perpendicular to the slope.

- 1 Lock the brakes on the C-arm system.

The red indicator appears (brakes locked).

- 2 Set the steering lever to lock position for transverse travel to the right or left, so that the wheels are braked in transverse position.

3.2.8 Transporting the monitor trolley

Preparing the monitor trolley


Prior to transport, the rolled-up power cable and the connection cable for the C-arm system should be placed on the handles of the monitor trolley.

CAUTION

Increased tipping effect.

Risk of crushing!

- ◆ Only products/components approved by Siemens may be installed.


- ◆ Hang the rolled-up cables on the handles of the monitor trolley.

Avoid transporting the monitor trolley on inclined surfaces $> 5^\circ$ or it may slip or tip over.

Moving the monitor trolley


- ◆ Release the central brake on the front of the monitor trolley (see arrow) to start moving it.

When the monitor trolley is moved, care must be taken that it does not collide with, for example, catheters or anesthesia tubes.

Locking the monitor trolley

- ◆ To lock the monitor trolley, push the center brake down with your foot until it stops/engages.

3 System Description


CAUTION

Brake failure.

Risk of crushing!

- ◆ Before beginning the examination, perform the daily function and safety checks.

Setting the direction of the monitor trolley

- ◆ Push the center brake upwards with your foot until it stops/engages.

This locks the front right wheel into a fixed position when moving the monitor trolley, and thereby sets the direction of movement.

The monitor trolley then can be pushed in the direction of the set front wheel.

Monitor trolley park position

CAUTION

Insufficient braking force.

Risk of crushing!

- ◆ When parking the monitor trolley on an inclined surface, face the front toward the incline and lock the brakes.
- ◆ Make sure the monitor trolley is not parked on a wet floor or other places where the adhesion of the wheels is greatly reduced.


On inclined surfaces with an angle of $> 5^\circ$:

- 1 Make sure that the brakes are locked.
- 2 Make sure that the front of the monitor trolley faces forward in the direction of the incline.

Park position of the monitors with height-adjustable monitor column (optional)

The monitor trolley holds two rotatable monitors. Using the corresponding option, these can be height-adjusted and folded up for transport.

(1) Height-adjustable monitors

(2) Toggle switch for height-adjustment of the monitors

The height of the monitors can be adjusted with the toggle switch.

1 Press the arrow symbol for Down on the toggle switch.

Motors move the monitors down; they can only be lowered to a specific end position.


CAUTION

Motorized movements.

Risk of crushing!

- ◆ Do not stand in the travel range of the monitor or place objects on the trolley housing.

2 Prior to transport fold the monitors together to prevent colliding with obstructions.

3 System Description

CAUTION

Folding the monitors.

Risk of crushing!

- ◆ Be careful when folding in the monitors - your hands could get caught.

3.3 Emergency Situations

Keep a replacement unit ready for critical applications in an emergency. Refer to (→ Page 19 *System failure*).

3.3.1 Function of the Emergency Stop button

A red EMERGENCY STOP button is located on the C-arm system control unit and the (optional) remote control.

The EMERGENCY STOP button stops the following processes immediately in a hazardous situation:

- motorized C-arm movements
- active radiation

When the EMERGENCY STOP button is activated (pressed), a red lock symbol flashes on the active emergency stop button. The lock symbol lights constantly on the other emergency stop components as long as the emergency stop button is active in at least one location.

If the EMERGENCY STOP button has been activated (pressed), radiation can still be released manually.

- ◆ Unlock the EMERGENCY STOP button by turning it to the right only when the hazardous situation has been eliminated.

3.3.2 Behavior in case of power failure

In the case of a power failure, the unit is switched off. In order to prevent data/image loss, an uninterruptible power supply (UPS) which carries out a controlled shutdown of the imaging system is installed in the monitor trolley. The images from the last examination are saved if the auto save function is switched on (**Auto save** must be on for interventional applications).


The following are not saved:

- Collimator position
- Set dose
- Image orientation
- Position memory
- Zoom

If there is one, a hospital emergency power supply cuts in if the line power supply fails. However, this usually occurs with interruption and can take some time.

The system can be restarted as soon as the mains power supply is available or the emergency power unit is functioning.

Restarting

- 1 Press the ON button on the monitor trolley to restart the system.
- 2 Wait until the unit has fully rebooted.
- 3 Pay attention to error messages.

When the power-up process has been completed (after approx. 2 minutes), the Cios Fusion is fully operational again.

- 4 Reselect the last patient in order to continue his or her treatment.
- 5 Perform the procedure, dose and collimator settings required in order to continue the application.

The positioning of the C-arm remains intact when the system is switched off. Nevertheless, please check the current position and correct it if necessary.

3.3.3 Restarting the Cios Fusion

The Cios Fusion must be restarted following certain error messages.

The following settings are reset to their default values by a restart:

- Collimator position
- Set dose
- Image orientation
- Position memory
- Zoom

- 1 Switch the system off by pressing the **OFF** button on the monitor trolley.
- 2 Wait until the computer has shut down.

- 3 Press the **ON** button on the monitor trolley to restart the unit.
- 4 Pay attention to the error messages displayed.
- 5 Reselect the last patient in order to continue his or her treatment.
- 6 Perform the procedure, dose and collimator settings required in order to continue the application.

The positioning of the C-arm remains intact when the system is switched off. Nevertheless, please check the current position and correct it if necessary.

CAUTION

Restarting the imaging system is not possible.

Treating patients no longer possible!

- ◆ Restart the system as presented in the following reset procedure.


3.3.4 Restart via reset

If the computer fails to shut down or no longer reacts to inputs, a system reset may be necessary.

The following information will be lost following a reset:

- Collimator position
- Set dose
- Image orientation
- Position memory
- Zoom

Resetting

- 1 Use a pin or a pointed object to actuate the reset button on the monitor trolley (the system will not automatically boot up afterwards).
- 2 After the system switches off, press the ON button to switch the unit on. Wait until the system has rebooted.
- 3 Reselect the last patient in order to continue his or her treatment.
- 4 Perform the procedure, dose and collimator settings required in order to continue the application.

The positioning of the C-arm remains intact when the system is switched off. Nevertheless, please check the current position and correct it if necessary.

3.3.5 Procedure for cardiopulmonary resuscitation (CPR)

If cardiopulmonary resuscitation (CPR) becomes necessary, the Cios Fusion must be removed from the OR table to ensure free access to the patient.

- 1 Move the C-arm into a position where it can be removed from the table without causing a collision.
- 2 Release the brakes and move the C-arm to the required position either via motor drive or manually. Reapply the brake to prevent unwanted movements of the C-arm.
- 3 If necessary, remove the remote control unit from the patient table and fasten it to the monitor trolley.
- 4 Then release the brake lever and move the Cios Fusion out of the treatment area.
- 5 At the same time, watch out for the cables leading to the monitor trolley, the remote control unit and the footswitch.
- 6 Release the brake of the monitor trolley, if necessary, and remove it from the treatment area.

3 System Description

4 Examination

4.1	Safety information	97
4.1.1	Prior to the examination	97
	Examination settings	97
	Audible warnings	97
	Software failure	97
	Storage capacity	97
	Reference images from previous examinations	97
4.1.2	During the examination	98
	Patient position	98
	Image orientation	98
	Data conformity	98
	Detector temperature	98
4.2	Preparing the system	99
4.2.1	Removing/inserting the scattered radiation grid	99
	Removing the scattered radiation grid	99
	Inserting the scattered radiation grid	100
4.2.2	Protection against contamination and the penetration of fluids	100
4.2.3	Positioning the C-arm	100
	Using the laser light localizer	100
	Aligning the C-arm	101
4.2.4	Unlocking/locking radiation release	101
	Unlocking radiation release	101
	Locking radiation release	102
4.2.5	Selecting the application	102
	Setting the application and application group	103
4.3	Patient registration	103
4.3.1	Emergency registration	105
	Emergency registration on the control panel	105
	Emergency registration on the monitor	105
	Emergency registration with the hand switch or footswitch	105
	Provisional patient data	106
4.3.2	Registering a new patient	106
	Opening patient registration	106
	Data Entry Dialog window	106
	Entering data	108
	Completing registration	109
4.3.3	Registering previous patients	109
	Opening the Patient list	109
	Checking/correcting patient data	110
	Completing registration	110
	Resuming a commenced study	111
4.3.4	Registering patients from the worklist	111
	Updating the worklist	111
	Opening the worklist	112
	Checking/correcting patient data	113
	Completing registration	113
	Deleting a patient entry	113

4 Examination

4.4	Defining the examination settings	114
4.4.1	Changing the fluoroscopy parameters	114
	Selecting the operating mode	114
	Setting collimators	116
	Selecting the dose rate level	117
	Switching on metal correction	118
	Selecting the motion/noise filter	118
	Selecting the input format (Detector zoom/ magnification)	119
	Setting the pulse rate	119
	Setting the X-ray parameters manually	120
4.4.2	Setting the image parameters	121
	Selecting the image orientation	121
	Standard image orientation	122
	Rotating/flipping images	122
4.5	Acquiring images	123
4.5.1	Releasing radiation	123
	Displaying fluoroscopy images	124
	Switching on enlarged preview image (blow-up)	125
	Confirming a warning signal	125
	Temperature monitoring	126
	Radiation information	126
4.5.2	Saving and displaying images	126
	Storing images (during radiation)	127
	Storing images (after radiation)	127
	Saving and reviewing a scene	127
	Changing the image display	128
	Recording video	129
4.5.3	Using reference images	129
	Transferring Images	129
	Scrolling through images	130
	Holding the reference image	130
	Loading images from an external video source	130
4.6	Ending the examination	131
4.6.1	Closing the patient	131
4.6.2	Examining the next patient	132
	Another patient registered	132
4.7	Performing special examinations	133
4.7.1	Subtraction angiography (Sub)	133
	Progression	133
	Performing subtraction angiography	133
4.7.2	Roadmap (Road)	134
	Fill image available	134
	Generating a new fill image	134
	Using the fill image from subtraction angiography	135
	Manually selecting the fill image	135
	Positioning the catheter (phase C)	136
4.7.3	Using graphics overlays	137
	Live graph overlay menu	137
	Drawing lines	137
	Deleting points	137
	Postprocessing lines	138
	Deleting everything	138
4.7.4	Pediatric use	138

4.1 Safety information

4.1.1 Prior to the examination

Examination settings

Before starting surgery, please make sure that all the set parameters as well as examination settings in all operation modes are correct.

Audible warnings

Audible warnings over the loudspeakers indicate possible hazards, such as crushing hazards, longer radiation release or high dose rate.

CAUTION

An audible signal sounds when below the safety standard.

Risk of crushing!

- ◆ Check for proper loudspeaker function before beginning the examination. Note when booting the system whether an audible signal is heard.

Software failure

CAUTION

Software crash.

Working with the imaging system is no longer possible.

- ◆ Restart the system.

Storage capacity

The hard disc fill level is displayed in the status bar of the left monitor. A warning signal is displayed before the final capacity of the hard disk is reached.

This icon, for instance, indicates that about 25 percent of the disk space has been used.

This icon, for instance, indicates that the available disk space has reached a critical level.

Please make sure that there is sufficient storage capacity before you start the examination. In addition, please observe the relevant system messages.

Reference images from previous examinations

CAUTION

Risk that the network connection may not be available continuously.

Reference images not available during intervention.

- ◆ Before beginning the intervention, the user should load the required images from the network archive.

4.1.2 During the examination

CAUTION

C-arm movement during surgery.

Risk of injury!

- ◆ Due to the maneuverability of the C-arm, collisions with other surgical devices are possible if the device is not used as intended and in accordance with the instructions in this manual.

CAUTION

Instruments falling into the gap between the control unit and its housing.

Risk of short circuit, electric shock, and/or fire!

- ◆ Make sure no instruments, conductive material, or fluids penetrate into the openings/joints between the control unit and its housing.

Patient position

Prior to the release of radiation, it is necessary to check whether the patient is positioned correctly.

Image orientation

Ensure correct image orientation (accurate to side) on the monitor/on the hardcopy.

Data conformity

Correct orientation of image and patient data must be checked before storage.

Before terminating an examination and beginning work on the next patient, the patient data should be verified.

The registered patient should be deselected at the end of the examination.

(→ Page 131 *Closing the patient*)

Detector temperature**CAUTION**

Detector temperature is not within recommended range.

Risk of limited image quality in the resulting X-ray images. Risk of X-ray radiation exposure without any diagnostic value!

- ◆ Regularly check the temperature icon on the display.

Exceeding the critical temperature threshold of 54 °C is indicated by a system message on the left monitor. In this case it is recommended to acknowledge the message and shut down the system immediately to allow it to cool.

4.2 Preparing the system

⚠ CAUTION

Brake failure.

Crushing hazard from automatic movements!

- ◆ Before beginning the examination, perform the daily function and safety checks.

4.2.1 Removing/inserting the scattered radiation grid

Removing the scattered radiation grid reduces the radiation dose (as well as image quality), e.g., for a pediatric examination.

⚠ CAUTION

Changing the grid.

Crushing hazard between the detector housing parts.

- ◆ Note the hazardous zones/locations between the detector housing parts.


⚠ CAUTION

Incorrect grid status: inserted, removed.

Radiation unsuitable for imaging results in unnecessary radiation exposure!

- ◆ Select the application according to your clinical procedure.
- ◆ Make sure that the current grid status meets your image requirements. The imaging system indicates whether a grid has been inserted or removed.
- ◆ If necessary, insert or remove the grid manually.

Removing the scattered radiation grid

- 1 Loosen the knurled screw (1) until you can pull the grid out.
- 2 Carefully pull the scattered radiation grid from its holder (2).
- 3 Put the scattered radiation grid in a safe place.

- 4 Check that the icon for "Grid inserted" has a line through it.

This indicator appears on the left monitor and the control panels in the exposure parameter area.

Inserting the scattered radiation grid

Once the pediatric examination is completed, you must reinstall the scattered radiation grid in order to ensure optimum image quality during standard examinations.

- 1 Install the scattered radiation grid in the reverse order of its removal: Push in the grid and tighten the knurled screw.
- 2 Check that the icon for "Grid inserted" does not have a line through it.

This indicator appears on the left monitor and the control panels in the exposure parameter area.

4.2.2 Protection against contamination and the penetration of fluids

If a significant amount of fluid is expected during an examination, there is a risk of fluids penetrating into the system. It is recommended to cover the relevant areas appropriately.

The C-arm can be covered completely or partially with a sterile disposable cover to protect it against contamination (→ Page 296 *Sterile cover on the C-arm*).

CAUTION

Fluid penetration and contamination.

Damage to the flat detector

- ◆ Make sure that the detector is always covered by a disposable cover during the examination.

4.2.3 Positioning the C-arm

Using the laser light localizer

Depending on the configuration of your Cios Fusion you can use the laser light localizer as a positioning aid.

Switching on

- ◆ Press this button on the control panel.

Depending on the system configuration, all laser light localizers on the C-arm system switch on.

The target crosshairs are displayed on the left monitor as long as the laser light localizer is switched on (depending on the configuration settings).

Switching off

- ◆ Press the activated button on the control panel.

All activated laser light localizers are switched off automatically after 5 minutes.

Aligning the C-arm

CAUTION

Collision with other systems.

Risk of injury!

- ◆ Due to the maneuverability of the C-arm, collisions with other systems are possible if the device is not used as intended and in accordance with the instructions in this manual.

- 1 Align the C-arm system to the table and patient.
- 2 Release the relevant brake and move the C-arm to the required exposure position (→ Page 67 *Move C-arm*).
- 3 Lock all brakes again.

4.2.4 Unlocking/locking radiation release

Your system has a lock-out option for radiation release to prevent unintentional release of radiation, e.g., by accidentally pressing the footswitch.

The radiation lock can be optionally secured with a password so only authorized persons can set the exposure parameters and release radiation.

You can set and change the password for locking radiation release in the configuration (→ Page 219 *Password for radiation release*).


Unlocking radiation release

- 1 Press this button on the **PREPARATION** task card.

Without password protection on the radiation lock, exposure parameters can be set and radiation released.

With radiation lock password protection, the password entry dialog is displayed.

- 2 Enter the correct password (default setting: 0000).
Exposure parameter setting and radiation release are released.
– or –

- Press this button if you want to keep the lock in place.
The dialog box closes without changes.

If you mistype your password, you can delete your entry with the **Del** button.

Locking radiation release

- 1 Press this button on the **PREPARATION** task card.
Without password protection on the radiation lock, exposure parameter setting and radiation release are locked.
With password protection, the password entry dialog is displayed.
- 2 Enter the correct password.
Exposure parameter setting and radiation release are locked.
– or –

- Press this button if you do not want to lock radiation release.
The dialog box closes without changes.

4.2.5 Selecting the application

The applications contain individual sets of examination settings tailored to various medical application areas (such as General, Ortho).

The available applications are divided into application groups, for example, according to hospital-specific requirements.

CAUTION

The image orientation is reset.

Risk of an incorrect decision during surgery!

- ◆ Be aware that the image orientation is reset when changing the application group or application. All previous settings are lost.

For further information about applications and application groups, as well as their configuration, see (+ Page 219 *Examination settings (PEX Editor)*).

Setting the application and application group

The application groups and applications are available in selection lists on the **PREPARATION** task card.

Application group	General/Std.	▼
Application	Standard Pat.	▼

- 1 Press the arrow button to open the **Application group** or **Application** selection lists.

- 2 Use the arrow buttons to scroll up and down in the selection lists.

- 3 Click the desired entry.

The exposure parameters including the operating mode are preset according to the selected application.

4.3 Patient registration

Before you can examine a patient with your system, you must register the patient.

Registration means that you give your system all the information about a patient that it requires for an examination.

Do not press the key combination **Ctrl + Alt + Del**.

CAUTION

Risk of mistaking patients.

The wrong patient is being treated.

- ◆ The user must ensure the correct patient is being treated before he or she is registered on the system.

CAUTION

Patient ID already exists.

Misinterpretation of data of the examined patient if the selection is based solely on the patient ID.

- ◆ Always use four-way identification [patient name, sex, date of birth, patient ID].
- ◆ Avoid mixing up patient data by using four-way identification. Selecting patient data based on the patient ID alone is insufficient. Ensure the proper display of parameters in the **Worklist** and in the **Patient list**.

Depending on how registrations are organized in your hospital and how much time you have for registration, you can choose between different patient registration procedures.

■ Emergency registration

If a patient is admitted who is in an extremely critical condition and must therefore be examined and treated immediately, select emergency registration. The patient is registered with provisional data as an emergency patient. This reduces the time before you can begin the examination to a minimum.

■ Registration for the examination

However, if you want to register a patient for an examination, you first enter the patient's data or select it from the database and then examine the patient.

■ Preregistration

If you want to prepare the system to examine a patient at a later point in time, then you can preregister the patient.

For example, in the morning you can enter the data of all the patients to be examined during the day. When you want to begin an examination, simply select the relevant data and edit them, if necessary. This saves time during the examination.

■ HIS/RIS query

If your system is connected to a HIS/RIS system (hospital and radiology information system), the patient can be registered for the examination through one of these systems.

CAUTION

HIS/RIS configuration not harmonized.

Communication not possible.

- ◆ Within the HIS/RIS do not use the "Only Patient ID" configuration because the patient is identified in the imaging system based on four pieces of data (patient name, ID, sex, date of birth).

Patient registration is also possible on the monitor trolley when the C-arm system is disconnected.

4.3.1 Emergency registration

An emergency registration is performed if a patient must be examined and treated immediately, without any time to waste entering the patient's data.


Emergency registration on the control panel

Emergency registration can be performed on the control panel at any time, regardless of the processing step you find yourself in.

PREPARATION task card

- 1 Change to the **PREPARATION** task card if necessary.

- 2 Press the "Emergency patient" button.

The **PROCESSING** task card is displayed on the control panel.

You can now start with the examination.

Keyboard mode

If you are in a processing step requiring text entry, keyboard mode is set on the control panel.

- ◆ Press this button in the symbol button bar.

The patient is registered as an emergency patient.

Keyboard mode ends.


Another patient registered

If you perform emergency registration while another patient is registered, a corresponding dialog box is displayed.

- 1 Select the desired options to conclude the examination (→ Page 131 *Closing the patient*).
- 2 Press **OK** to close the currently registered patient.

Emergency registration on the monitor

When you are editing patient data the **Data Entry Dialog** window is displayed on the monitor. You can start emergency registration in that window with a click of the mouse.

- ◆ Click this button in the **Data Entry Dialog** window.

The patient is registered as an emergency patient.

The dialog window on the monitor closes.

Emergency registration with the hand switch or footswitch

Outside a running examination, you can perform emergency registration by activating the hand switch or footswitch.

- 1 Make sure that there is not a patient registered for an examination.
- 2 On the hand switch, press the radiation release key **once**.

– or –

On the footswitch, press the radiation release pedal **once**.

The patient is registered as an emergency patient.

Provisional patient data

During emergency registration, missing patient data are supplemented automatically. An emergency patient is registered with the following data:

Data type	Entry
Last name	Emergency
First name	Emergency if the patient ID has already been entered sequential number if the patient ID has not already been entered
Date of birth	Today's date
Sex	Unknown
Patient ID	E_sequential number

Make sure to correct the provisional patient data when the opportunity becomes available (→ Page 146 *Correcting data*).

4.3.2 Registering a new patient

If a patient has never been examined in your hospital or practice before, no data about this patient will be stored on your system. Therefore all the required data for this patient must be entered before an examination.

Opening patient registration

Open patient registration on the control panel in the **PREPARATION** task card.

1 Change to the **PREPARATION** task card if necessary.

2 Press this button.

The control panel switches to keyboard mode.

The **Data Entry Dialog** window opens on the left monitor.

Data Entry Dialog window

The mandatory fields for registration and preregistration are highlighted in bold. Which fields are displayed and which are mandatory can be configured (→ Page 209 *Patient registration*).

PATIENT area Name, Patient ID, Age, and Sex of the patient are always mandatory entry fields. This information uniquely identifies the patient in your databases.

Special notes:

■ **Patient ID**

You can have a unique patient ID generated as long as no internal rules have been set for the format.

■ **Date of Birth**

Four digits must be entered for the year of birth.

If the date of birth is unknown, the patient's age is used to calculate a date of birth from today's date.

■ **Age**

If the date of birth has been entered, the patient's age is calculated from today's date and entered.

STUDY area In this area, data from the HIS/RIS system (optional), if connected, are entered automatically. The data are used especially to uniquely identify the study created for the examination.

INSTITUTION area This area contains information regarding the examining facility and personnel. The information can be helpful if the examination images are passed on to a different organization for reporting.

Special notes:

■ **Institute Name**

The institute name entered during configuration is entered here as the default.

Application area This area shows the applications available for selection, divided into different application groups (→ Page 102 *Selecting the application*).

Entering data

The cursor is in the input field for the patient's name.

Making entries

- 1 Enter text and numbers using the keyboard in the control panel.

The characters entered are shown on the monitor and in the text field on the control panel.

- 2 Use the mouse to select entries from the selection list.

You can manage the content of certain selection lists (drop-down combo boxes) in the following way:

- To add an entry to the list, mark the text entered in the text field and press **Shift + ↑** (arrow up) on the keyboard.
- To remove an entry from the list, select it from the list (text is marked in the text field) and press **Shift + ↓** (arrow down) on the keyboard.

Changing fields

Tab

- ◆ Use the **Tab** button (hold down the **Shift** button to go backwards) to switch to the required field.

Generate a patient ID

- ◆ Press this button in the symbol button bar on the control panel.


– or –

Click the icon in the **Data Entry Dialog** window.

A new patient ID is generated based on the current date and time.

Correcting the entry

- 1 On the control panel press the arrow button of the text field to move the cursor character-by-character.

Back space

- 2 Press this button to delete the character in front of the cursor.

– or –

Shift

Mark the incorrect characters with the arrow buttons while holding the **Shift** button, and overwrite with the keyboard.

Deleting all entries

- ◆ Press this button in the symbol button bar on the control panel.
All entries in the **Data Entry Dialog** window are deleted.


Completing registration

Once you have at a minimum completed the mandatory fields in the **Data Entry Dialog** window, you can register the patient for the ensuing examination or preregister him or her for examination later on.

Examining the patient

- ◆ Press this button in the symbol button bar on the control panel.
– or –

Click this button in the **Data Entry Dialog** window.

The **Data Entry Dialog** window closes.

The **PROCESSING** task card is displayed on the control panel.


You can start examining the patient.

Preregistering patients

During preregistration the patient and entered data are stored locally. You can then access the patient data again when you start the examination (→ Page 111 *Registering patients from the worklist*).

- ◆ Press this button in the symbol button bar on the control panel.
– or –

Click this button in the **Data Entry Dialog** window.

The patient is included in the **Worklist**.

The **Data Entry Dialog** window closes.

The **PREPARATION** task card is displayed on the control panel.

4.3.3 Registering previous patients

A patient who has already been examined in your hospital or practice is known to your system.

For the current patient registration you apply the saved information in the **Data Entry Dialog** window.

Opening the Patient list

You open the **Patient list** on the control panel in the **PREPARATION** task card.

- 1 Change to the **PREPARATION** task card if necessary.

- 2 Press this button.

The control panel switches to keyboard mode.

The **Patient list** opens on the left monitor.

★ Patient Name	Patient ID	Date of Birth	Sex	Study Description	Physician	Flag
Emergency, 61	E_61	16-01-2013 (0)	O	Standard		-----
Emergency, 66	E_66	16-01-2013 (0)	O	Standard		-----C--
Kunz, Andrea	160113_104...	01-04-1945 (67)	F	Standard		-----C--
Meier, Hans	160113_103...	23-12-1976 (36)	M	Standard		-----C--
Muster, Bernhard	160113_104...	12-05-2011 (1)	M	Standard		-----C--
Schmidt, Manuel	160113_104...	12-12-1955 (57)	M	Standard		-----C--
Zander, Stefanie	160113_104...	09-05-2000 (12)	F	Standard		-----C--

Patient not found If the patient you want is not displayed, a database filter may be on or the patient has been deleted.

- 1 Set the filter criteria accordingly or switch off the filter. See (→ Page 145 *Filtering data*)
- 2 If necessary, import the patient data from the archive. See (→ Page 148 *Importing data*)

Checking/correcting patient data

Opening data

- 1 Use the arrow buttons or mouse to select the required patient.

- 2 Press this button in the symbol button bar on the control panel.
The data of the selected patient are transferred to the **Data Entry Dialog** window.

Editing data

- ◆ Check that the information in the **Data Entry Dialog** window is correct and add any missing data (→ Page 108 *Entering data*).

Wrong patient

If you have selected the wrong patient:

- ◆ Press this button in the symbol button bar on the control panel.
All changes in the **Data Entry Dialog** window are discarded.
The **Patient list** is displayed again.

Completing registration

Examining the patient

- ◆ Press this button on the control panel or click the button on the screen.
The **Data Entry Dialog** window closes.
The **PROCESSING** task card is displayed on the control panel.
You can start examining the patient.


Preregistering patients

- ◆ Press this button on the control panel or click the button on the screen.
The patient is included in the **Worklist**.
The **Data Entry Dialog** window closes.
The **PREPARATION** task card is displayed on the control panel.

Resuming a commenced study

If you cancelled a study with **Close**, you can continue this study later. Just select the study, load the data into the **Data Entry Dialog** window, and then start the examination. All newly generated images are saved in the existing study as a new series. Please note that this is only possible if you have already acquired images in the commenced study. An empty study cannot be used for registration.

- 1 Open the **Patient list**.
- 2 Select the corresponding study.
- 3 Transfer the data to the **Data Entry Dialog** window.

- 4 Register the patient for examination.

The **PROCESSING** task card is displayed on the control panel.

You may resume your study.

4.3.4 Registering patients from the worklist

The **Worklist** contains patients that were preregistered manually in your system and automatically via the HIS/RIS connection.

Updating the worklist

The **Worklist** is updated via the HIS/RIS system at regular intervals, and when the **Worklist** is opened. You can also initiate the update manually, for example, when you want to continue working without HIS/RIS connection.

- ◆ Press this button on the **PREPARATION** task card on the control panel.

Manage Tool is displayed on the left monitor.

Worklist complete

- ◆ Click this button in the **Patient** subtask card.

The **Manage Tool** closes.

The **Worklist** opens on the left monitor and is updated with the examinations planned for your system.

Filtering the worklist

If you are looking for a specific patient you can enter his name or you can query the HIS/RIS for specific entries only.

4 Examination

- 1 Click this button in the **Patient** subtask card.

The **Get Worklist Dialog** window opens.

The **Get Worklist Dialog** window contains the following fields:

- PATIENT**
 - Last Name: [Text Field]
 - First Name: [Text Field]
 - Patient ID: [Text Field]
- INSTITUTION**
 - Ref. Physician: [Dropdown Menu]
 - Ref. Physician: [Dropdown Menu]
- STUDY**
 - Accession No.: [Text Field]
 - Request ID: [Text Field]
- Modality**
 - ☐ CR ☐ Other
 - ☐ DX
 - ☐ XA
 - ☐ RF
- ☐ Not only today's worklist
- Buttons:** Get Worklist, Cancel

- 2 Enter the data you know as search criteria.

Possible placeholders: * stands for one or more random characters, ? stands for one random character.

Modality: The modality planned for the patient.

Not only today's worklist: Query for all scheduled appointments.

Get Worklist

- 3 Click this button.

The **Manage Tool** closes.

The **Worklist** opens on the left monitor and updates. In addition to the patients preregistered manually, only those HIS/RIS entries that meet the search criteria are displayed.

Opening the worklist

You open the **Worklist** on the control panel in the **PREPARATION** task card, unless it is already displayed due to a manual update.

PREPARATION

- 1 Change to the **PREPARATION** task card if necessary.

Worklist

- 2 Press this button.

The control panel switches to keyboard mode.

The **Worklist** opens on the left monitor.

Patient Name	Patient ID	Date of Birth	Sex	Study Description	Procedure Step	Locatl...
Emergency, 65	E_65	16-01-2013 (0)	O			
Meier, Hans	160113_103444	11-11-1967 (45)	M			
Meier, Hans	160113_103359	11-11-1986 (26)	F			
Emergency, 64	E_64	16-01-2013 (0)	O			

Checking/correcting patient data

Opening data

- 1 Use the arrow buttons or mouse to select the required patient in the **Worklist**.

- 2 Press this button in the symbol button bar on the control panel.
The data of the selected patient are transferred to the **Data Entry Dialog** window.

Editing data

- 1 Check that the information in the **Data Entry Dialog** window is correct and add any missing data (→ Page 108 *Entering data*).

- 2 Press this button on the control panel or click the button on the screen to confirm changes.

The **Worklist** is displayed again.

Correct Patient

For reasons of data consistency, some information in HIS/RIS entries cannot be changed, such as the personal data of the patient.

Wrong patient

If you have selected the wrong patient:

- ◆ Press this button in the symbol button bar on the control panel.
All changes in the **Data Entry Dialog** window are discarded.
The **Worklist** is displayed again.

Completing registration

- ◆ Use the arrow buttons or mouse to select the required patient in the **Worklist**, if necessary.

Examining the patient

- ◆ Press this button on the control panel to examine the patient.
The **Worklist** closes.

The **PROCESSING** task card is displayed on the control panel.

You can start examining the patient.


– or –

Double-click the patient with the mouse.

Deleting a patient entry

As soon as the examination of a patient in the **Worklist** ends, the data saved in the study is moved automatically to the **Patient list**.

If an examination for a patient is cancelled, you can manually delete the entry so that the **Worklist** remains clear.

1 Use the arrow buttons or mouse to select the required patient in the **Worklist**, if necessary.

2 Press this button on the control panel to remove the patient in question from the **Worklist**.

4.4 Defining the examination settings

The standard exposure parameters are already set by default for the examination, depending on the **Application** you have selected (→ Page 102 *Selecting the application*).

If these settings are appropriate for your examination, you can start the image acquisition immediately on the Cios Fusion (→ Page 123 *Acquiring images*).

The **Application** currently selected together with the patient data are displayed in the left monitor on the upper left.

Meier, Hans
23-12-1976
Standard

Make sure that the patient's name on the monitor matches the patient to be examined.

4.4.1 Changing the fluoroscopy parameters

On the control units, you set the exposure parameters in the left area of the control panel. Settings may be changed prior to each radiation release.

CAUTION

Incorrect exposure parameters.

Risk of unnecessary radiation exposure due to incorrect exposure parameters.

- ◆ Ensure the exposure parameters are correct. Pay attention to the exposure index (EXI).

Selecting the operating mode

The operating mode currently set is indicated by the button highlighted white. The current operating mode is deselected when switching to another mode.

The operating modes can also be selected via the corresponding button on the multifunctional footswitch (optional), if used. You can then release radiation again right away. (→ Page 291 *Multifunctional footswitch (optional)*)

Fluoroscopy

Fluoroscopy operating mode enables real-time imaging which, for example, allows viewing of moving objects.

- ◆ Press this button.

Fluoroscopy mode is selected.

The button is highlighted white.

Upon switching on the Cios Fusion, the operating mode configured as the standard application, usually Fluoroscopy, is selected automatically by default (→ Page 219 *Examination settings (PEX Editor)*).

Single image

Single image operating mode enables an electronic instant image with the best image quality. It is recommended for final exposures.

On activation of single image a short radiation pulse is released.

Please note that a Fluoro exposure with > 2fps has to be made before a single exposure in order to set the exposure parameters (kV/mA) and thereby achieve good image quality.

When acquiring single images, remember to keep the radiation release button pressed (on the hand switch or footswitch) until the radiation indicator goes out.

- ◆ Press this button.

Single image operating mode is selected.

The button is highlighted white.

After complete image acquisition, radiation is automatically switched off, even if the radiation release button remains pressed.

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The subtraction technique enables an isolated display of the vascular system after injection of the contrast agent by means of background subtraction (→ Page 133 *Subtraction angiography (Sub)*).

- ◆ Press this button.

Subtraction operating mode is selected.

The button is highlighted white.

Roadmap

The Roadmap technique enables the user to position a catheter precisely in a blood vessel under fluoroscopy. (→ Page 134 *Roadmap (Road)*)

- ◆ Press this button.

The **Roadmap** mode is selected.

The button is highlighted white.

Please note that a Fluoro exposure with > 2fps has to be made before exposures using subtraction or roadmap technique in order to set the exposure parameters (kV/mA) and thereby achieve good image quality.

Setting collimators

The Cios Fusion uses the following collimators to collimate the X-ray beam:

■ Rectangular collimator

The rectangular collimator provides radiation protection for the patient and all persons participating in the examination.

■ Slot collimator

The slot collimator is used primarily for collimation of the extremities.

⚠ CAUTION

Wrong collimator setting.

Dose not as required!

- ◆ The collimators need to be set so that at least one edge of the leaves is visible in the fluoroscopy image. This must be checked daily.

Digital collimator function: The areas outside the collimated areas are automatically blackened if so configured for the selected application.
(→ Page 225 *Editing applications*)

Collimation menu

- ◆ Press this button.

The **Collimation** menu opens.

The current collimation setting is displayed on the left monitor and in the image area of the control panel.

When you release radiation, the collimators are in the position shown in the image.

When you open/close the rectangular collimator or move the slot collimator without radiation, you can see the position of the collimator on the LIH image displayed with lines superimposed.

Opening/closing the rectangular collimator

When switching on the Cios Fusion, the rectangular collimator automatically opens to full format. Smaller collimation produces less scatter radiation and therefore better image contrast.

- ◆ Press one of these buttons.

The rectangular collimator opens/closes.

Opening/closing the slot collimator


Collimation enhances image contrast and reduces scatter radiation. Direct radiation that passes the soft tissue laterally is reduced to such an extent, that differences in brightness do not disturb when images are viewed on the monitor.

The position of the left and right slot collimator blades can be changed jointly (symmetric) and individually (asymmetric).

- 1 Press one of these buttons.

The slot collimator opens/closes symmetrically.

- 2 Press one of these buttons for the left or right slot collimator side.
The corresponding blades extend/retract on one side.

Rotating the slot collimator

By rotating the slot collimator, the collimated field can be quickly oriented to the direction of the anatomy under examination (e.g. the extremities).

- ◆ Press one of these buttons.
The slot collimator rotates to the left/right.

Resetting the collimators

- ◆ Press this button.
The collimators are reset to full format.

Selecting the dose rate level

You can choose between three dose levels (deviating from the current user program).

The dose level currently set is displayed at the **Dose** button of the control panels and at the left monitor.

When selecting different dose levels, please note that additional specific parameter settings, such as k factor and pulse rate, change automatically.

Dose menu

- ◆ Press this button.
The **Dose** menu appears.
The current setting is highlighted white.


Selecting the dose level

- ◆ Select the required dose with the buttons.
Dose levels of low, medium, and high may be set.

If the maximum skin dose is exceeded at the highest dose level, the corresponding button is called "High level". The maximum duration of a fluoro scene is limited to 30 s. The button and the display on the monitor are marked with a warning symbol.

Switching on metal correction

Metal correction is designed for examinations with metal parts in the measurement range (e.g., larger, non-standard implants). When activated, the contrast and brightness settings are adjusted accordingly and the automatic dose rate control is improved.

- ◆ Press this button.
The button is highlighted white.

Selecting the motion/noise filter

By using noise filters with motion detection, you take into account the influence of movements when filtering noise. Image integration depends on the preset k factor (depending on the operating mode up to max. $k = 16$), which is adapted according to the selected motion/noise filter.

k factor: a number of k exposures are integrated into one image; the k factor can be selected between $k = 1$ and $k = 16$ and can be assigned to an examination set and stored (IEC 60601-2-7:1998 29.1.103 d).

Motion noise filter menu

- ◆ Press this button.
The **Motion Noise filter** menu is displayed.
The current settings are highlighted white.

Motion detection

Motion detection is used to establish whether the maximum number of images (for still objects) or a lower number of images (for moving objects) should be used for image integration.

- ◆ Press one of these buttons.
Off: Image integration is independent of movements in the image with preconfigured k-factor.
On: Image integration takes place only if movements are detected in the image.

Cardiac filter

The **Cardiac** button is used for examinations requiring particularly high temporal resolution. It reduces the k factor to a minimum to reduce blurring. It may however also increase noise.

- ◆ Press this button.
No image integration occurs. In this way rapidly occurring processes (e.g. beating heart) can be best shown.

Averaging live image

By fine adjusting the averaging, you take into account the motion frequency of the object to be examined. An integration factor is applied to the k value that results for the selected motion detection.

- ◆ Press a button in the **Noise filter live image** area.
Low: A low integration factor is selected (for recording fast-moving objects; default: 0.7).
Medium: No change in the k value (integration factor = 1).

High: A higher integration factor is selected (for very slow movements; default: 1.4).

Noise filter LIH

A separate setting for image integration is available for the LIH (Last Image Hold) image. The values to be set and the functionality is the same as for **Noise filter live image**.

- ◆ Press a button in the **Noise filter LIH** area.

Possible settings: **Low, Medium, High**.

Selecting the input format (Detector zoom/magnification)

Three different zoom levels are available. This sets the input format for the subsequent image acquisition. The input format is the size of the detector area that is exposed and then fills the screen when displayed on the monitor.

Zoom level	Mag 0	Mag 1	Mag 2
Input format (20 cm x 20 cm detector)	20 cm x 20 cm (full format)	15 cm x 15 cm	10 cm x 10 cm
Input format (30 cm x 30 cm detector)	30 cm x 30 cm (full format)	20 cm x 20 cm	15 cm x 15 cm

Magnification menu

- ◆ Press this button.

The **Magnification** menu opens.

The current input format level is displayed on the left monitor and in the image area of the control panel.

Select the zoom level

- ◆ Select the required zoom level with the keys.

The input format display is updated accordingly.

Setting the pulse rate

Depending on the configuration level, you can select from among up to 9 different pulse rates during the examination. The respective set value is displayed on the control panel using the **Pulse rate** button.

Pulse rate menu

- ◆ Press this button.

The **Pulse rate** menu opens.

The current setting is highlighted white.

Selecting the pulse rate

- ◆ Select the pulse rate required with the corresponding button.

The pulse rate changes accordingly.

Setting the X-ray parameters manually

You can change the default X-ray parameters manually with the +/- buttons for kV/mA or kV/mAs (for single image operating mode).

Automatic dose rate control: Using automatic dose rate control (ADR), the kV/mA values or kV/mAs values are regulated so that the mean value of the image gray values is kept constant largely independently of the object thickness and position. This ensures optimal image quality for on-screen evaluation.

The current automatic dose rate control status can be read on the control panel using the **Tech lock** button.

- (1) ADR activated (kV/mA are automatically regulated)
- (2) ADR locked (kV/mA remain constant for the current values)

The system possess an automatic metal correction which is designed for examinations with metal parts in the measurement range. The metal correction parameters are included in the **Application** settings. For example, the contrast and brightness settings are adjusted accordingly and the automatic dose rate control is improved.

If the image quality is not sufficient, please use the additional **Metal correction** or **Tech lock** functionality.

Activating Tech lock Stop

When metallic objects (e.g. intermedullary nails) are introduced into the beam path or when examining objects of varying density (e.g. hip prosthesis) under fluoroscopy, it is recommended that you set the kV just established with the **Tech lock** button at the start of fluoroscopy.

The automatic dose rate control (ADR) must remain activated during interventional procedures.

- 1 Press this button.
The **Tech lock** menu opens.

- 2 Press this button.
The stop function is activated and the **Stop** button is highlighted white.
Automatic dose control is disabled.
The +/- buttons for kV/mA or kV/mAs are enabled.

Do not switch between Fluoroscopy and Single Image operating modes while **Tech lock Stop** is active, because with **Tech lock Stop** the image brightness cannot be adjusted to the different dose levels in the two operating modes.

Changing the X-ray parameters manually

Activating **Tech lock Stop** enables you to manually set kV/mA or kV/mAs. The set values are displayed on the control panel.

- (1) Display of kV/mA (mAs) and pulse width in the **Tech lock** menu
 (2) Display of kV/mA (mAs) in the exposure parameter area

- ◆ Press the +/- buttons.

The values for kV/mA or kV/mAs increase or decrease.

Pressing the button for an extended period continually increases or decreases the corresponding X-ray parameter.

Once the upper or lower limit of the setting range is reached, an audible tone is emitted. Pressing the button again has the same effect.

The mA values assigned to the kV values result from the fluoroscopy curves (→ Page 249 *Curves and diagrams*).

4.4.2 Setting the image parameters

On the control units, you set the image parameters in the **PROCESSING** task card in the right area of the control panel. Changing the settings impacts the image currently displayed (if applicable) and any images acquired thereafter.

Selecting the image orientation

It may be necessary to rotate and/or flip the image so that it appears in the required orientation on the monitor.

Object display on the monitor depends upon the C-arm system position relative to the patient.

CAUTION

Inverted image orientation.

Risk of an incorrect decision during surgery.


- ◆ The operator must apply due diligence; identifying each side is recommended.

CAUTION

Inverted image orientation.

Incorrect intervention decision!

◆ It is recommended to mark the patient position using lead letters.

Standard image orientation

(1) Right side of patient

(2) Left side of patient

Rotating/flipping images

The current image rotation and -flip settings are displayed in the image area of the control panel.

U: 0° not flipped

The rotation angle ranges from 0° to 359°.

Rotation Flip menu

◆ Press this button.

The **Rotate Flip** menu appears.

Selecting the image rotation

1 Use these buttons to set the angle of rotation in 90° increments.

– or –

Press one of these buttons.

With each short press of the button, image rotation increases 1° in the corresponding direction.

Pressing the button for an extended period continually increases or decreases the angle of rotation.

2 Press this button to reset the angle of rotation to 0°.

Selecting image flip

◆ Press one of these buttons.

Vertical/horizontal image flip is activated.

4.5 Acquiring images

All further steps of the examination are performed at the C-arm system. For detailed information on operating the C-arm, hand switch, and footswitch, see (→ Page 49 *System Description*).

Meier, Hans
23-12-1976
Standard

The registered patient's last name, first name and date of birth appear in the upper control area of the left monitor.

◆ Prior to releasing radiation, make sure that the registered patient and actual patient are the same.


CAUTION

Inverted image orientation.

Risk of an incorrect decision during surgery.

◆ Prior to releasing radiation, the user must ensure the correct patient position with respect to the X-ray beam.

4.5.1 Releasing radiation

◆ Release radiation with the hand switch or footswitch.

The current fluoro image is displayed on the left (live) monitor.

Please keep the radiation release switch pressed until radiation has been released and a usable image is displayed.

Complete image integration as a function of the set k factor (up to max. k = 16, depending on the operating mode) is ensured even for very short exposures ("toe tapping").

k factor: a number of k exposures are integrated into one image; the k factor can be selected between $k = 2$ and $k = 16$ and can be assigned to an examination set and stored (IEC 60601-2-7:1998 29.1.103 d).

CAUTION

Defective imaging system hardware.

Risk of X-ray radiation exposure without a saved image!

- ◆ If no image appears, immediately stop the radiation and switch off the tube assembly.

CAUTION


Because of a technical error the acquisition is not stopped or there is an unintended release of radiation.

Risk of unnecessary radiation exposure!

- ◆ Press the **OFF** button on the monitor trolley, even if the **EMERGENCY STOP** button has already been pressed. Radiation terminates immediately, even if the sound of the rotating anode can still be heard.

Displaying fluoroscopy images

In the image area the fluoroscopic images are displayed during and after exposure. Additional information may be shown as image text, depending on the configuration (→ Page 214 *Image display and print properties*).

- (1) Information on kV and mA
(2) Window value and LUT information

Live images As soon as the exposure (i.e. radiation) starts, the current fluoroscopic images (live images) are displayed on the left monitor and the control panels.

Last Image Hold (LIH) The last image hold (LIH) is displayed as soon as the exposure has ended.


Switching on enlarged preview image (blow-up)

You can always enlarge (blow up) the displayed image using the control panel.

1 Tap the image area.

The displayed image is enlarged across the entire task card range.

Buttons to change the image display appear on the right side.

(1) Blow-up area

(2) Image display functions

2 Press the Close button to reset the control panels to normal display.

Confirming a warning signal

After 5 minutes of uninterrupted fluoroscopy time, a warning signal sounds. However, this can be reset on the control panels. After 10 minutes of uninterrupted fluoroscopy time, radiation is terminated. However, it can be released again.

These settings are country-specific and can be changed by Siemens Customer Service in accordance with the applicable regulations.

◆ Press the reset button (alarm icon flashes).

The audible warning signal is deactivated.

The fluoro timer is reset.

4 Examination

Temperature monitoring

If the system is used intensively over an extended period, the X-ray tube and single tank can become very hot.

Single-tank temperature

The current thermal capacity of the single tank is shown by an 11-segment bar display on the left monitor:

The length of the bar indicates the relative temperature load of the X-ray system.

X-ray tube temperature

The X-ray tube temperature status is displayed on the left monitor as follows:

If the temperature is normal, a green tube symbol appears on the color monitor.

If the tube is heated, an orange tube symbol appears on the color monitor.

If the tube temperature is critical, a red tube symbol with an exclamation point appears on the color monitor.

To prevent overheating, the pulse rate for a very hot tube is limited to 10 images per second. Working with the reduced frame rate is nevertheless possible without any further restrictions.

Radiation information

The radiation data are shown on the control panels and on the left monitor. The information is updated with every release of radiation.

00:00:11
hh:mm:ss

Here, the entire fluoroscopic time since the start of the examination of a patient is displayed.

55
μGym²

The cumulative dose area product for the current patient is displayed if the optional dose measurement chamber is installed.

As an alternative: Display of the air kerma value and cumulative air kerma value. These data appear on the left monitor.

The reference location for determining the air kerma values with the optional dose measurement chamber is 30 cm in front of the detector input. The reason for this convention is that in typical applications the object to be examined is located approximately 30 cm in front of the flat detector.

Depending on country-specific regulations, you can have Siemens Customer Service change the display so that it indicates the air kerma and cumulative air kerma value instead of the dose area product.

4.5.2 Saving and displaying images

In Fluoroscopy, Subtraction, and Roadmap operating modes, images can be saved manually during the examination. Images that you want to print out later have to be saved during the examination. Images printed directly during the examination are saved automatically.

If the application used is configured such that all images are saved automatically, manual saving is not necessary and therefore not possible.

Saving images/scenes with the footswitch is also possible when the pedal allocation is configured accordingly.

Images/scenes can always be saved via the corresponding button on the multifunctional footswitch (optional). See (→ Page 291 *Multifunctional footswitch (optional)*).

Storing images (during radiation)

- ◆ Press this button on the control panel.

– or –

Press this button on the hand switch.

The current image (or the subtraction image in Sub/Roadmap mode) is stored in the local database. All images saved in the course of an examination are stored as a series.

The saved fluoro image is displayed on the right (reference) monitor. This does not occur if a second video source is displayed on the right (reference) monitor (refer to (→ Page 130 *Loading images from an external video source*)).

Storing images (after radiation)

- ◆ Press this button on the control panel.

– or –

Press this key on the hand switch (press and hold < 2 seconds).

The displayed image (LIH) is stored in the local database.

The LIH is displayed on the right (reference) monitor.

Holding the key for > 2 seconds: Saves the scene last recorded (LSH).

Saving and reviewing a scene

In the fluoroscopy, subtraction and Roadmap operating modes, scenes can be saved and reviewed (Last Scene Hold) during the examination.

Saving the scene

- ◆ Press this button on the control panel.

– or –

Press this key on the hand switch (press and hold > 2 seconds).

The current scene (LSH) will be stored in the local database. In case of a restricted review range, images are saved from the starting point to the end point.

The LSH is displayed on the right (reference) monitor.

4 Examination

Holding the key for < 2 seconds: Saves the image last recorded (LIH).


Reviewing a scene

- 1 Press this button.
The **Cine replay** menu appears.

- 2 Press the activated Pause button to start (replay speed matches the storage rate).
During replay the Pause button changes to a Start button.


- 3 Replay at half/normal/double speed.

- 4 Previous/next single image.

- 5 Play previous/next available scene.

- 6 Establish starting point and end point for replay.
The corresponding image number is displayed in the respective button.

- 7 Switches automatic cine replay on/off.

If **Autoreplay** is configured in the application used, cine replay is started automatically at the end of the exposure.

Changing the image display

There are tools available on the control panels for subsequent optimization of the image display. For detailed information see (→ Page 152 *2D image processing*).

- 1 Edge enhancement and spatial noise filter (4 levels each).

- 2 Brightness/Contrast, image inversion.

- 3 Zoom/pan image/set shutter.

- 4 Rotate/flip an image.

CAUTION

Processing images by changing the image display.

Loss of image quality by artifacts.

♦ Use care when changing the image display settings.

Recording video

Parallel to the display on the left monitor, you can save scenes to DVD as video. The resulting DVD is suitable for use with a DVD player.

Starting recording


- 1 Insert an empty DVD into the DVD drive.

Supported formats: DVD-R, DVD+R

After inserting a DVD when a recorder license is available, the monitor queries whether the DVD should be used for recording.

- 2 Change to the **EXTERNAL** task card.


- 3 Press this button.

Ending the recording

- 1 Change to the **EXTERNAL** task card if necessary.

- 2 Press this button.
3 After the write process, remove the DVD from the drive.

4.5.3 Using reference images

In addition to the live images on the left monitor, selected images, e.g. images used for comparisons, can be displayed on the right monitor.

Transferring Images

Automatic

During saving, an image or scene is transferred automatically to the right monitor:

- 1 Save the image manually with the **Store** button (during/after radiation)

The displayed image/LIH is saved and transferred.

– or –

Save the image automatically ("**Automatic save**" is active)

The LIH is saved and transferred after radiation ends.

- 2 Save the scene manually with the **Store cine** button

The LSH is saved and transferred.

– or –

Save the scene automatically ("**Automatic save**" is active)

The LSH is saved and transferred after radiation ends.

Manual

Using the **Transfer A to B** function you can transfer an image to the right monitor during or after radiation, without saving it simultaneously.

- ◆ Press this button.

The displayed image is transferred.

Please remember that the transferred image has not been saved, and will be lost permanently if overwritten by another image.

Scrolling through images

If the "**Automatic save**" function is active or you have saved several times, images in addition to the displayed image are available on the right monitor.

- 1 Scroll to previous image/to previous scene.

- 2 Scroll to next image/to next scene.

Holding the reference image

The **Hold reference** function allows you to select a particular image to permanently display as the reference image on the right monitor.

- 1 Press this button.

The display on the right monitor is "frozen".

The transfer of additional images and scrolling to other images is disabled on the right monitor.

- 2 Press the button again to cancel **Hold reference**.

Hold reference is automatically turned on for the first image transferred to the right monitor using the **Transfer A to B** button.

Hold reference is cancelled automatically when switching to image overview display (→ Page 153 *Displaying the image overview*). Hold reference is reestablished once the image overview display is closed.


Loading images from an external video source

If suitable video sources are connected to your system (such as endoscope or microscope), you can display the recorded images on the right monitor for reference.

Selecting the video source


1 Change to the **EXTERNAL** task card.

2 Press this button.

The video inputs of your system are checked according to the connecting video sources.

3 Keep the button pressed until the required video source is reached.

Displaying images

◆ Press this button.

The image recorded from the video source is displayed on the left monitor if radiation is not being released.

The image is saved to the local database as "CAP image" under the current study of the registered patient.

4.6 Ending the examination

After you have completed your exposures, the last acquired image is displayed on the left monitor. Now you finish the examination of the current patient.

If you want to examine the next patient immediately afterwards, you can register that patient straight away. In this case the examination of the current patient is automatically finished.


4.6.1 Closing the patient

Closing the patient ends the current examination. The patient data and images saved in a study in the local database are moved to the **Patient list**.

1 Press this button on the bottom right of the control panel.

The **Close** menu appears.

2 Select the desired options.

Print rad. summary: Prints the radiation report on the default printer.

Export rad. summary: Sends the radiation report to the default address for export.

Send MPPS: Sends the documentation to the HIS/RIS or a network address established for this purpose.

3 Press this button.

The examination is concluded. The patient data and images are cleared from the monitors.

– or –

Press this button to continue the examination with the registered patient.

4.6.2 Examining the next patient

1 Change to the **PREPARATION** task card if necessary.

2 Register the next patient for examination (→ Page 103 *Patient registration*).

Another patient registered

If another patient is still registered the **New patient** menu is displayed.

A blue dialog box titled 'New patient' with a warning icon. The text inside says: 'Another patient is still registered. Close current patient with below checked options and continue with new patient?'. There are three checkboxes, each with a white checkmark: 'Print rad. summary', 'Export rad. summary', and 'Send MPPS'. At the bottom are two buttons: 'OK' and 'Cancel'.

1 Select the desired options to conclude the examination (→ Page 131 *Closing the patient*).

2 Press **OK** to close the currently registered patient.

The new patient is registered for examination.

4.7 Performing special examinations

The following optional operating modes require special procedures deviating from the standard examination:

- Subtraction angiography (Sub)
- Roadmap (Road)
- Using live graphic overlays
- Pediatric applications

After incorrect operation, it may be necessary for the user to repeat certain procedures, e.g. administering of contrast agent.

4.7.1 Subtraction angiography (Sub)

During the examination, images without contrast agent (mask) are continuously subtracted from images with contrast agent and displayed on the monitor. Depending on the contrast agent flow, they display the relevant vascular segment without superimposition in real time.

Progression

Subtraction angiography is divided into three phases:

- Phase A
Time until the mask is completed (permanently defined)
- Phase B1
Time from the "inject" display on the monitor until the contrast agent has reached the examination region
- Phase B2
Time of the actual exposure of the examination region


Performing subtraction angiography

- ✓ The patient to be examined has been registered.
- 1 Make the necessary settings for the examination. (→ Page 114 *Defining the examination settings*)
- 2 Press this button.

Subtraction operating mode is selected.


Phase A

- ◆ Release radiation with the hand switch or with the assigned footswitch.

During the generation of the mask the native image is displayed on both monitors.

Phase B 1 Keep the exposure release button pressed.

2 Inject the contrast agent as soon as the syringe symbol appears on the screen and control panel.

The corresponding subtraction image is displayed on the control panel and left monitor.

In the native image on the right monitor you can see the continuous filling of the blood vessel with contrast agent.

Keep the radiation release button pressed until the vessel is filled with contrast agent.

Depending on the setting of the **Landmark** function, a certain percentage of anatomical background is displayed in the subtraction images (→ Page 227 *Parameter for fluoroscopy and exposure*).

Fill image

♦ Let go of the radiation release button.

After radiation has been switched off the maximum fill image is calculated and displayed on the left monitor.

The vessels are highlighted in the native image on the right monitor.

The mask (phase A), subtraction series (phase B) and the fill image are saved automatically.

4.7.2 Roadmap (Road)

In the first step of the Roadmap examination, the fill image is created as a mask from a normal subtraction (phase A and phase B).

In the second step, the display of the vessel into which the catheter is to be positioned is superimposed by current fluoroscopic images (phase C).

Fill image available

You can use the existing fill image if you have already performed a subtraction angiography during the current examination (→ Page 135 *Using the fill image from subtraction angiography*).

The availability of the fill image is shown on the **Road** button and on the footswitch icon on the monitor. In that case, the first step of the Roadmap examination is not required.

A previously created fill image can also be used if you have switched to a different operating mode (such as single image).

Generating a new fill image

If there is no suitable fill image (e.g. from a previous subtraction angiography), you have to generate a new one.


✓ The patient to be examined has been registered.

1 Make the necessary settings for the examination. (→ Page 114 *Defining the examination settings*)

- 2 Press this button.
The **Roadmap** mode is selected.

Phase A

- ◆ Release radiation with the hand switch or with the assigned footswitch.
During the generation of the mask the native image is displayed on both monitors.

Phase B

- 1 Keep the exposure release button pressed.
- 2 Inject the contrast agent as soon as the syringe symbol appears on the screen.
In the native image on the right monitor you can see the continuous filling of the blood vessel with contrast agent.
The corresponding subtraction image is displayed on the left monitor.

Keep the radiation release button pressed until the vessel is filled with contrast agent.

Fill image

- 1 Let go of the radiation release button.
Radiation is stopped.
The fill image is calculated and displayed.
- 2 Start positioning the catheter (→ Page 136 *Positioning the catheter (phase C)*).

Using the fill image from subtraction angiography

- ✓ Subtraction angiography (**Subtraction** operating mode) was performed at some point during the current examination.

- 1 Press this button.
The **Roadmap** mode is selected.
- 2 Start positioning the catheter (→ Page 136 *Positioning the catheter (phase C)*).

If you press the **Sub** button, the existing fill image is discarded. You will then have to regenerate the fill image.

Manually selecting the fill image

- ✓ A subtraction data set (subtraction or Roadmap operating mode) was generated in the current examination.


Fill image with another mask

To calculate a new fill image you have to determine a suitable subtraction mask that can be used.

- 1 Press this button.

The **Select mask** menu opens.

- 2 Move the slider to the left/right to select another mask image.

Pressing the arrow repeatedly takes you incrementally from image to image.
The fill image is recalculated automatically.


Fill image from image selection

To calculate a new fill image you have to determine a suitable selection of subtraction images that can be evaluated.

- 1 Press this button.

The **Fill image** menu opens.

- 2 Display the subtraction series images, and use the buttons to set the first/last image you select as the new start point/end point.

The **Fill image** menu displays the images numbers for the start point and end point of the limited series range.

A new fill image is calculated and displayed based on the image selection.

- 3 Press this button to apply the fill image.

The **Fill image** menu closes.

Selecting a different fill image

If multiple fill images have been calculated you can select the most suitable image for phase C.

- 1 Press this button.

The **Fill image** menu opens.

- 2 If necessary, repeatedly press one of the buttons to select the fill image required.

Positioning the catheter (phase C)

If the vessels are easily discernible, position the catheter in the fill image under fluoroscopy.

- 1 Release radiation with the hand switch or with the assigned footswitch.

- 2 Position the guidewire or the catheter under fluoroscopic control.

The right monitor shows the fluoro image (native).

The left monitor shows the subtracted Roadmap image with catheter.

You can repeat fluoroscopy as often as needed while you insert the catheter.

4.7.3 Using graphics overlays

The lines drawn in an ongoing examination can be transferred to the subsequent live images. Such graphic overlays are used, for example, in orthopedics, traumatology as well as in subtraction and roadmap examinations. For example, lines can be drawn into fill images to support the treatment of abdominal aortic aneurysms.

WARNING

An annotated image with graphic overlay lines is displayed. Following system and/or table or patient movements and/or image modifications, subsequent images may not match the overlay line displayed due to projective changes.

Incorrect interpretation of the display and resulting incorrect decision during operation!

- ◆ The graphic overlay line is displayed until the user actively turns it off. This is done independently of any changes (for example, movements, operating mode changes, image modifications). The user must ensure that the displayed overlay line matches the current, underlying image (for example, display of anatomical structures via the label function and/or use of contrast agent).

Live graph overlay menu

- 1 Display the required image.
- 2 Press this button.

The Live graph overlay menu is displayed.

Drawing lines

You draw the lines as a polygon where each change in direction creates a vertex in the image.

- 1 Press this button to set the starting point and vertices at the respective mouse pointer positions.

- 2 Press this button at the required end point.

The new line is shown as a dotted line in the image.

If a graphic overlay was already used in the previous image, the corresponding lines are displayed as continuous. The old overlay is replaced by the newly drawn in lines in the following images.

Deleting points

If a point has been placed incorrectly, it can be deleted immediately after being drawn as long as the end point has not been set.

- ◆ Press this button.

The last point set is deleted from the image.

Pressing the button again deletes the previous point.

You can continue drawing the line.

Postprocessing lines

Postprocessing via the **Live graph overlay** menu is possible only for newly drawn lines (dotted). A previously used overlay (solid lines) can be post-processed using the **Measure Annotate** menu see (→ Page 166 *Evaluating images*).

- 1 Move the mouse pointer to the line or a point to select the given graphic object.

- 2 Press this button.

Now you can move (arrow buttons) or delete (**Delete selected** button) the selected point or the complete line.

Deleting everything

- ◆ Press this button.

The newly drawn graphic overlay (dotted lines) is deleted.

4.7.4 Pediatric use

The Cios Fusion may be used for pediatric purposes.

We suggest using the available special applications for pediatric use in the application database delivered with the system. There are individual pediatric applications for each application group of the default application database. These default pediatric applications may be used for basis of the examination for all pediatric patients (in general for children up to the age of 14 years).

It is recommended to remove the scattered radiation grid for all pediatric patients up to the age of 7 years as well as all medical applications that are restricted to extremities.

For pediatric patients older than 7 years a scattered radiation grid is recommended otherwise image quality may be compromised.

If depending on patient stature the kV increases above 100 kV it is suggested to switch to the corresponding adult program (otherwise image quality may be compromised).

Depending on the medical requirements of your respective examination, consider other dose saving measures. For example, reduce the frame rate and radiation time as much as possible.

Following these suggestions will lead to a minimized dose for the patient.

For an orientation of dose savings please see the following table with a comparison between pediatric and standard protocols:

Object: 10 cm PMMA directly located on the FD cover (without scattered radiation grid for pediatric programs)

Application (Group)	Fluoro		Sub	
	Standard - Dose rate (mGy/min)	Pediatrics - Dose rate (mGy/min)	Standard - Dose rate (mGy/min)	Pediatrics - Dose rate (mGy/min)
Card	1.4	0.2	1.6	0.4
Gastro	1.4	0.2	1.5	0.2
General	0.8	0.2	1.1	0.2
Ortho	0.8	0.2	1.1	0.2
Uro	0.8	0.2	1.1	0.2
Vascular	0.8	0.2	1.0	0.2

5 Postprocessing

5.1	Patient data	143
5.1.1	When to use the Patient list	143
5.1.2	External data	143
5.1.3	Searching for and selecting patient data	143
	Flag	144
	Opening the Patient list	144
	Sorting data	145
	Filtering data	145
	Selecting data	146
5.1.4	Correcting data	146
	Opening the correction dialog	146
	Data Entry Dialog window	146
	Editing data	147
5.1.5	Deleting data	148
	Deleting data manually	148
	Protecting data from deletion	148
5.1.6	Importing data	148
	Import dialog	149
	Opening the import dialog	149
	Data from removable devices	149
	Data in the network	150
5.2	2D image processing	152
5.2.1	Loading and displaying images	152
	Loading images	153
	Displaying the image overview	153
	Changing the display mode	154
	Scrolling	155
	Displaying full screen image	155
	Reviewing a scene	156
5.2.2	Selecting images	157
	Single selection	157
	Multiselection	157
5.2.3	Deleting individual images/series	159
5.2.4	Changing the image display	160
	Selecting images for processing	160
	Windowing images	160
	Inverting images	161
	Using image filters	162
	Rotating/flipping images	163
	Zooming and panning images	163
	Setting digital shutters	164
	Resetting the image display	165
5.2.5	Evaluating images	166
	Selecting images for evaluation	166
	Displaying the evaluation functions	166
	General work steps	167
	Calibration	168
	Measuring the distance	169
	Measuring an angle	170
	Annotating images	170
	Setting the side identification	171

5 Postprocessing

5.2.6	Evaluation in enhanced mouse mode	171
	Activating enhanced mouse mode	172
	Turning off enhanced mouse mode	172
	General work steps	172
	Calibration	173
	Measuring the distance	173
	Measuring an angle	173
	Annotating images	174
5.2.7	Saving images	175
	Saving a single image	175
	Saving multiple images	175
5.3	Subtraction processing	176
5.3.1	Postprocessing functions	176
5.3.2	Loading and displaying images	176
	Loading images	176
	Changing the monitor layout	177
	Displaying images	177
5.3.3	Changing the image display	177
	Windowing images	177
	Including anatomical background (Landmark)	178
	Making the image and mask coincide exactly (Pixel shift)	178
5.3.4	Evaluating subtraction images	180
	Stenosis quantification	180

5.1 Patient data

All acquisitions (images, series) prepared as part of an examination are saved in a study on the hard drive of the imaging system (local database). You can access these data quickly and easily using the **Patient list**.

5.1.1 When to use the Patient list


- To examine a patient who has already been examined once before with your system and whose data are still saved in the local database
(→ Page 109 *Registering previous patients*)
- To view the images of a patient from earlier examinations in order to compare them with current results
- To postprocess or add comments to images after an examination
- To correct incorrect information on a patient stored in your system
- To print a patient's images for documentation purposes
- To archive patient data and images or to send them to another location in your hospital via the network

5.1.2 External data

To view or process data from an external storage medium or other network node, they first have to be imported.

5.1.3 Searching for and selecting patient data

The studies of patients examined previously that are saved to the local database are displayed in the **Patient list**.

- (1) Table heading
- (2) Study entries (List of patients and their studies)
- (3) Preview images of the selected study

Flag

An abbreviation in the **Flag** column indicates the study status:

Flag	Meaning
A/a	All/some images were sent to the archive and the receipt was confirmed (as long as the archive works with DICOM Storage Commitment)
C	Study is closed (no further acquisition possible)
D/d	All/some images were transferred to or from CD/DVD
P/p	All/some images were printed
R/r	All/some images read in via network
S/s	All/some images were sent and the receipt was confirmed (if the recipient works with DICOM Storage Commitment)
#	Study is protected against deletion
*	Study was received from HIS/RIS

Opening the Patient list

You open the **Patient list** on the control panel in the **PREPARATION** task card.

 PREPARATION

- 1 Change to the **PREPARATION** task card if necessary.

 Previous patients


- 2 Press this button.

The control panel switches to keyboard mode.

The **Patient list** opens on the left monitor.

Sorting data

You can sort list entries according to the displayed columns.

 Patient Name

- 1 Click the table header for the respective column.

The list is sorted in alphanumerically ascending order based on this column.

Studies with the same entry in the sorting column are further sorted by the date and time of the examination.

- 2 Click the table header of the same column again.

The sorting sequence is reversed.

Filtering data

It is possible to only display studies that meet specific criteria by using database filters. The **Flags** of the individual studies are evaluated in this case (→ Page 143 *Searching for and selecting patient data*).

Activating a filter

 Filter

- 1 Press this button in the symbol button bar on the control panel.

The **Database Filter** window appears on the left monitor.

- 2 Activate the required check boxes.

show: only display studies with the selected flag

hide: hide all studies with the selected flag

- 3 Confirm with **OK** to apply the filter.

5 Postprocessing

You can specifically delimit the list display by combining multiple filter criteria. Example: **Printed show** together with **Archived hide** shows the studies that have been printed but not yet archived.

Deactivating the filter

- 1 Press this button in the symbol button bar on the control panel.

The **Database Filter** window appears on the left monitor.

- 2 Click this button.

All data (unfiltered) are displayed in the **Patient list**.

Selecting data


Searching by the first letter

- 1 Using the keyboard in the control panel, press the first letter of the required patient.

The first name with a suitable patient name is selected in the **Patient list**.

- 2 Press this same letter button again, to select the next list entry for a patient with this first letter.

Direct selection

- ◆ Use the arrow buttons on the control panel to select the required study.

– or –

Click the list entry with the mouse.

5.1.4 Correcting data

Every now and then it is necessary to correct patient data for an examination. This is the case particularly with emergency patients that are saved to the local database with provisional data.

Patient data in studies you have already exported or sent via the network cannot be changed in the local database.

Opening the correction dialog

Patient and examination data are corrected in the **Data Entry Dialog** window.

- 1 Select the study entry to be corrected in the **Patient list**.

- 2 Press this button in the symbol button bar on the control panel.

The data of the selected patient are transferred to the **Data Entry Dialog** window.

Data Entry Dialog window

Patient and examination data of the selected study saved in the local database are entered in the **Data Entry Dialog** window.

Editing data

Which fields are displayed and therefore can be corrected depends on the configuration settings (→ Page 209 *Patient registration*).

Changing fields

Tab

- ◆ Use the **Tab** button (hold down the **Shift** button to go backwards) or click the mouse to switch to the required field.

Correcting the entry

- 1 On the control panel press the arrow button of the text field to move the cursor character-by-character.

Back space

- 2 Press this button to delete the character in front of the cursor.
– or –

Shift

Mark the incorrect characters with the arrow buttons while holding the **Shift** button, and overwrite with the keyboard.

Saving changes

- ◆ Press this button on the control panel or click the button on the screen to confirm changes.

The **Data Entry Dialog** window closes.

The **Patient list** is updated based on the changes.

5.1.5 Deleting data

You should delete archived images you no longer need from the local database on a regular basis. This ensures sufficient storage capacity on your system.

You can delete data as follows:

■ Automatic deletion

Depending on the configuration settings, images are deleted once they reach a specific status (→ Page 218 *Deleting images automatically*)

■ Manual deletion

You can individually select studies and delete regardless of data status

Deleting data manually

1 Select the entry to be deleted in the **Patient list**.

2 Make sure that the selected study has been archived correctly.

3 Press this button in the symbol button bar on the control panel.

The selected study and all its images are deleted.

During image postprocessing, you can delete individual images or series of a study using the image overview (→ Page 159 *Deleting individual images/series*).

Protecting data from deletion

You can protect images on the local database against unintentional manual deletion and automatic deletion.

Setting deletion protection


1 In the **Patient list** select the entry you want to protect against deletion.

2 Press this button in the symbol button bar on the control panel.

Removing deletion protection

1 Select the deletion-protected entry in the **Patient list**.


2 Press this button in the symbol button bar on the control panel.

5.1.6 Importing data

If the required patient data and images are on removable devices or other network nodes, import them to the local database. The data are then available in the **Patient list**.

Import dialog

- (1) Data source selection
- (2) Display thumbnails (CD/DVD and USB storage media only)
- (3) Disk space required on the hard drive (CD/DVD and USB storage media only)
- (4) Button bar
 - Import:** Starts data import
 - Status:** Displays the data transfer status
 - Close:** Closes the dialog window

Opening the import dialog

- ◆ Click this button on the left monitor (lower left).
The **Import from...** dialog window is displayed.

Data from removable devices

Image data on a data medium in DICOM format can be read in as long as an appropriate drive is installed and configured.

Displaying saved data

- 1 Have the required data medium on hand.
 Insert the CD/DVD into the drive.
 Connect the USB memory to the USB socket.
- 2 Select the applicable data medium as the data source in the **Import from...** dialog window.
 The studies saved on the data medium are listed.
 If activated, thumbnails are displayed for the currently selected study.
 Example: Displaying the content of a CD/DVD

5 Postprocessing

Patient Name	Patient ID	Social Security ...	Military Rank	Date of Birth	Sex	Study Description
201, Emergency	E_201			18-10-2016 (...)	O	
200, Emergency	E_200			18-10-2016 (...)	O	Standard
199, Emergency	E_199			18-10-2016 (...)	O	General
198, Emergency	E_198			18-10-2016 (...)	O	General
197, Emergency	E_197			18-10-2016 (...)	O	standard
197, Emergency	E_197			18-10-2016 (...)	O	Standard
196, Emergency	E_196			18-10-2016 (...)	O	General
195, Emergency	E_195			18-10-2016 (...)	O	General
194, Emergency	E_194			18-10-2016 (...)	O	General
193, Emergency	E_193			18-10-2016 (...)	O	General
192, Emergency	E_192			18-10-2016 (...)	O	General
191, Emergency	E_191			18-10-2016 (...)	O	General
Patient, Road...	181016_1...			18-10-1972 (...)	O	General
paischent, nju	1234567			18-05-2016 (...)	F	General
Tech, Lock	181016_1...			18-10-1994 (...)	O	General
Store, Fluoro	181016_1...			18-10-1994 (...)	O	General
190, Emergency	E_190			18-10-2016 (...)	O	General
Fluoro, Modes	181016_0...			18-10-1961 (...)	M	General
fluoro, High s...	181016_0...			18-10-1994 (...)	F	General
Faktor, Kah	181016_0...			18-10-1983 (...)	F	Standard
189, Emergency	E_189			17-10-2016 (...)	O	Vascular
187, Emergency	E_187			17-10-2016 (...)	O	Standard
186, Emergency	E_186			17-10-2016 (...)	O	Standard
XA_AORTA	20120320...			10-10-1950 (...)	M	

Importing data 1 Click the required study with the mouse.

Press and hold the **Ctrl** button to select multiple entries.

Import

2 Click this button in the **Import from...** dialog window.

The data are read in.

Data in the network

You can access image data on other workstations or a long-term archive using "Query & Retrieve".

Opening "Query & Retrieve"

◆ In the **Import from...** dialog window select the **DICOM node** entry as the data source.

The **Query for retrieve** dialog window is displayed.

(1)

(2)

(3)

- The patients and studies found are listed in the results area.

5 Postprocessing

Patient Name	Patient ID	Birth Date	Sex	Study Name	Date and Time	Accession No.	Study Instance ID	Images	Series
Emergency	155893	12-11-1976	D		01-01-1900 12	Accession_35	15049105151		1
CR Modality	24886522	25-02-1923	D		03-12-2002 11	Accession_23	24861247221475		1
CR Modality	524752	14-06-1996	M		04-03-2005 09	abod	24814256542	40	1
Film	4745011	17-05-1911	F	Schädel	04-03-2006 02	abod	15111252		1
Service	276912		M	Schädel	04-03-2003 02	Accession_13	12453246547	10	1
Emergency 15	851855	14-08-1945	M	Schädel	04-03-2001 03	Accession_24	547541210	13	2

- (1) Details on the study currently selected in the results list
- (2) List of found patients and studies
- (3) Currently selected study

Displaying series

If you only require a selection of images from a specific study, then display the series contained in the study.

1 Select the corresponding study in the results list.

2 Click this button.

Series

All series of the selected study are listed in the results area.

Use this button to display the list of found patients and studies again.

Studies

Importing data

1 Select the required data (studies or series) in the results list with the mouse.

Press and hold the **Ctrl** button to select multiple entries.

2 Click this button in the **Import from...** dialog window.

Import

The data are read in.

5.2 2D image processing

As part of postprocessing, you view the results of an examination and evaluate them using 2D processing functions. You can also compare the image material from various examinations of a patient (such as the current examination and a previous examination).

5.2.1 Loading and displaying images

During an examination in progress you can only edit images of the registered patient. Outside of an examination, you can load the images of any patient for processing.

Loading images

In the **Patient list** you have access to the images of prior examinations stored on the local database. Select one or more studies of the required patient here (→ Page 143 *Searching for and selecting patient data*).

With an examination in progress, all saved images are already loaded. If needed, you can load additional images from prior examinations of the registered patient.

If the required images are no longer in the local database, then first import the data from the external source (→ Page 148 *Importing data*).

PREPARATION

Previous patients

Display
images

- 1 Change to the **PREPARATION** task card if necessary.
- 2 Open the **Patient list**.
- 3 Search for the patient, sorting and filtering the list as necessary.
- 4 Select the required study or studies.
- 5 Press this button in the symbol button bar on the control panel.

The selected images are loaded and displayed on the monitors.

The **PROCESSING** task card is displayed on the control panel.

– or –

Double-click the patient with the mouse.

Displaying the image overview

In the image overview display, the monitor layout changes specifically for 2D-postprocessing. You have direct access to all loaded images and you can select them for processing either individually or in any combination.

Activating the overview

If a patient is not registered the images are displayed in the image overview after loading. With an examination in progress, you may need to switch from the live image display to the image overview display.

Overview

- ◆ Press this button on the control panel.

– or –

On the left monitor click the 4x4 layout button (lower left).

The **Overview** menu is displayed on the control panel.

- (1) Selection buttons for Series/Images
- (2) Full image display
- (3) Multi-selection of Series/Images
- (4) Display mode (Stack/Stripe)

Monitor layout The following views are set on the monitors by default:

- Left monitor
 - 4x4 layout with all loaded images, sorted by exposure time
 - Stack display
 - 1st series selected (1st image segment is highlighted)
- Right monitor
 - 1x1 display of the (selected) 1st series

Switching off the overview

The image overview can only be switched off during an examination in progress.

- ◆ Press this button again.

The live image display is activated.

– or –

On the left monitor click the activated 4x4 layout button (lower left).

The system switches to live image display automatically during radiation release.

Changing the display mode

You can select from among the following displays of loaded images on the left monitor:

- Stack display: One series per image segment
- Image stripe display: One image from the selected series per image segment

Selecting stack display

Stack display is suitable as an overview of loaded studies and for direct comparison of images from several series.

- ◆ Press this button in the **Overview** menu.

Selecting image stripe display

- ◆ Press this button in the **Overview** menu.

In stack display, image segments that contain a scene are identified by an icon at the bottom right. Double-clicking the scene switches to image stripe display.

Scrolling

In the **Overview** menu you establish which image should be displayed on the right monitor.


Changing image/series

- 1 Select the suitable display mode:

Image stripe display if you want to show a different image in the current series.

Stack display if you want to display a different series.

- 2 Scroll in the image overview until the number of the required image or series is displayed.

- 3 Press the selection button for the required image or series.

The selected image is displayed on the right monitor.

– or –

In the left monitor click the applicable image segment with the mouse.

Next/Previous series

If you want to show images of another series in image stripe display, you can scroll through the series without switching to stack display.

- ◆ Press one of the buttons repeatedly if necessary.

Displaying full screen image

Basically, the image selected in the overview is displayed as a full image on the right monitor. However, if the **Hold reference** function is active the monitor display remains unchanged. In this case use the **Full image** function to display the required image in full size on the left monitor.

5 Postprocessing

Activating full image

- ◆ Press this button in the **Overview** menu.

Using the Next/Previous series buttons you can switch to other series in full image view.

Switching off full image

- ◆ Press this button again.

The image overview on the left monitor switches into the most recent display mode (stripe or stack).

Outside of an ongoing examination (no patient registered), you can turn full image view on and off using the activated 4x4 layout button (lower left of screen).

Reviewing a scene

After an examination, you can start movie mode to check the quality of the scenes you acquired.

A scene is usually reviewed on the right monitor. If **Hold reference** is active, the review takes place on the left monitor. The system switches automatically to 1x1 layout.

Cine replay menu

- 1 Show the required scene.

- 2 Press this button.

The **Cine replay** menu appears.

Starting the review


If the automatic review of scenes (**Auto replay**) is switched off, the review process is paused.

- 1 Press the activated Pause button to review the scene.

The replay speed varies with the number of saved images.


During review the Pause button changes to a Play button.

- 2 Press this button to change the review speed.

Replay at half/normal/double speed.

Limiting the playback range

- 1 Switch to the required image, for example using the slider or by stopping the review process at the appropriate location.

- 2 Press one of the buttons to set the starting point or end point for playback.

The corresponding image number is displayed in the respective button.

By pressing the starting point and end point buttons again you can reset the original playback range.

Scrolling image by image

- 1 Press this button to stop the review process.

- 2 Press the button to go backward/forward by a single image.

Different scene

- ◆ Press the button for the previous/next available scene.

5.2.2 Selecting images

You have to select images before you edit them. Single and multiple selections are possible.

- ◆ Select the suitable display mode in the **Overview** menu:

Image stripe display if you want to select multiple images in the current series.

Stack display if you want to select one or more series.

Single selection

- ◆ Press the selection button for the required image or series in the **Overview** menu.

– or –

In the left monitor click the applicable image segment with the mouse.

If the left monitor has been changed to full image mode, the image displayed there is selected.

Multiselection

Multiselection enables you to select images from one or more series individually. Images of a multiselection are marked on both monitors by a white dot.


- 1 Press this button in the **Overview** menu.

The menu extension for multiselection is shown.

The number of images currently selected is displayed.

5 Postprocessing

- Pressing the button again toggles to single selection.

Individual images/series

- Press the selection buttons for the required image or series in sequence.

— or —

Switch into keyboard mode, hold down the **Ctrl** button, and click the desired images/series on the left monitor.

Sequential images/series

- Press this button.

- Select the first required image/series.

- Select the last required image/series.
All images/series in between are selected as well.

— or —

Click the first desired image/series in the left monitor.

- Switch into keyboard mode, hold down the **Shift** button, and click the last desired image/series.

Supplementing the selection

- Switch the display mode if necessary.
- Select additional images/series using multiselection.

The newly selected images are added to the multiselection previously defined.

– or –

Press this button to add the image or series selected in the image overview to the multiple selection.

All images/series

◆ Press this button.

Image stripe display: All images in the series are selected.

Stack display: All loaded images are selected.

Deselecting images

◆ Deselect the selected images/series individually using the corresponding selection buttons in the image segments.

– or –

Press this button for the images or series to be deselected.

Press this button to deselect all images.

5.2.3 Deleting individual images/series

While the **Patient list** is used to delete complete studies (→ Page 148 *Deleting data*), you can use the image overview to select individual images or series for deletion.

1 In this case switch into stack display.

2 Switch into keyboard mode and press **Alt + d**.

A message box is displayed on the left monitor.

5 Postprocessing

3 Click the images/series to be deleted in the image overview.

The selected images are marked with an X.

4 Click **OK** to delete the selected images.

5.2.4 Changing the image display

In the loaded images you can adapt the window values (contrast and brightness), display enlarged sections, and center regions of interest in the middle of the image. You can also rotate, flip or invert images for certain diagnostic problems.

Image manipulation capabilities are limited for images imported from other modalities and images from external video sources (CAP images).

CAUTION

Incorrect or extreme image processing parameters.

Risk of limited image quality in the resulting X-ray images. Risk of incorrect decision during intervention or artifacts.

- ◆ Apply image processing parameters carefully, such as edge enhancement and windowing, and always check them carefully.
- ◆ Always check the image quality before reaching an interventional decision.

Selecting images for processing

Changes in image display affect all images of an individual series. The image functions cannot be applied to multiple series at the same time.

- ◆ Display an image of the required series (single selection).

Windowing images

Using the windowing settings, you establish a grayscale range in which the image is displayed. This enables you to highlight the region of interest and the relevant tissue type.

During windowing you optimize image brightness (shifting the windowing center along the grayscale) and image contrast (change of window width).

The set window values are shown at the bottom right in the image.

(1) WW: Window width – Contrast

(2) WC: Window center – Brightness

Brightness Contrast menu


- ◆ Press this button.

The **Brightness Contrast** menu is displayed.

Changing window values

- 1 Move the slider to the left/right to make infinitely variable adjustments in brightness (1) and contrast (2).


- 2 Click the left/right arrow to change the window values in small increments.

The controllers set numerical brightness values. Whether the image impression turns out brighter or darker depends on which display is selected (bone black or white).

Windowing using the mouse

- 1 Place the mouse pointer in the required image on the left monitor.

- 2 Move the mouse with the center mouse button pressed:

- (1) Brightness - (Center +)
- (2) Brightness + (Center -)
- (3) Contrast - (Width +)
- (4) Contrast + (Width -)

Fixed window value settings (LUTs)

Using window value settings, you can select three predefined settings for the window values using LUTs (look up tables). The values can be configured specifically for an application (→ Page 225 *Editing applications*).

- ◆ Press the button for the required preset.

Default

The **Default** setting assigns the window values from auto contrast.

- ◆ Press this button.

Brightness and contrast are reset to the default.

Inverting images

When an image is inverted, bright areas in grayscale images are displayed dark and dark areas are displayed bright. You use this function to toggle the bone display between white and black, for example.

5 Postprocessing

1 Open the **Brightness Contrast** menu.

2 Press this button.

Image inversion is switched on.

Using image filters

You can improve the appearance of blurred images using appropriate filters. The following filters can be used individually or in combination:

- Edge filter
- Spatial noise filter
- Edge-preserving noise reduction (optional)

Image enhancement menu

◆ Press this button.

The **Image enhancement** menu opens.

Selecting the edge filter

◆ Press one of these buttons.

The higher the filter setting, the more structures are enhanced and the less they are smoothed.

– or –

Press this button.

Edges are not enhanced.

Selecting the spatial noise filter

◆ Press one of these buttons.

The higher the filter setting, the more the image is smoothed, which will reduce its sharpness.

– or –

Press this button.

There will be no spatial noise filtering.

Selecting edge-preserving noise reduction (optional)

◆ Press one of these buttons (recommendation: start with **Low** and increase incrementally as needed).

The higher the filter setting, the more blurring is removed.

Large differences in gray values are not averaged so that the edges remain unaffected.

– or –

Press this button.

There is no edge-preserving noise reduction.

Rotating/flipping images

The rotate and flip functions allow you to easily compare images of series that were acquired in a different orientation.


Rotation Flip menu

◆ Press this button.

The **Rotate Flip** menu appears.

Selecting the image rotation


◆ Press one of these buttons.

With each short press of the button, image rotation increases 1° in the corresponding direction.

Pressing the button for an extended period continually increases or decreases the angle of rotation.

– or –

Use these buttons to set the angle of rotation in 90° increments.

Rotate using the mouse

- 1 Place the mouse pointer on the edge of the required image on the left monitor.
- 2 Drag the mouse pointer in the required direction of rotation.

No image rotation

◆ Press this button.

The angle of rotation is reset to 0°.

Selecting image flip

◆ Press one of these buttons.

Vertical/horizontal image flip is activated.

Zooming and panning images

You can zoom in on sections of images that you are particularly interested in (zooming). The zoomed image can be panned so that the relevant section is in the center of the image (panning).

Zoom Pan Shutter menu

- ◆ Press this button.

The **Zoom Pan Shutter** menu appears.

Zooming the image

- 1 Move the slider up/down for infinitely variable zoom in/zoom out of the image.

Upper limit: Zoom factor 4x

Lower limit: Zoom factor 1x

- 2 Click the up/down arrow to change the image size in small increments.

Panning the image

After zooming, if the relevant image area is outside the image segment, pan it accordingly.

- ◆ Press the arrow button for the direction in which the image content to be displayed is located.

Vertical arrow buttons: Pan up/down

Horizontal arrow buttons: Pan left/right

Diagonal arrow buttons: Combined horizontal and vertical pan

Image panning is possible only with a zoom factor > 1x.

Resetting Zoom/Pan

- ◆ Press this button.

The original image display is restored.

Setting digital shutters

You can hide irrelevant areas of the image. To do that, you place appropriate shutters over the image. Image areas covered by shutters appear black.

- 1 Open the **Zoom Pan Shutter** menu.

- 2 Press this button.

The menu extension for shutter definition is shown.

Adding shutters You can slide a total of 4 rectangular shutters into the image starting at the edges. First activate them individually or jointly using the appropriate buttons, and then fit them into place.

1 Press one of these buttons for the shutters you want to set.

Upper shutter

Right shutter

Left shutter

Lower shutter

2 Press this button to slide the activated shutters toward the center of the image.

3 Press this button to slide the activated shutters toward the edges of the image.

Hiding the shutters

◆ Press this button.

The full image area displayed again.

Resetting the image display

If adjusting the image parameters did not produce the desired results, you can reset the processing steps of the current session. All images are then returned to their original state.

- During an examination in progress: Status after image acquisition or after saving
- Outside an examination: Status after loading or after saving


After resetting, all changes in the image display are irrevocably lost.

1 Press this button.

A confirmation dialog box is displayed.

5 Postprocessing

2 Press **Yes** to reset the images.

– or –

Press **No** to close the dialog box without changes.

5.2.5 Evaluating images

Using graphical evaluation functions you can measure distances, lengths, and angles of relevant structures, as well as add annotations. The left monitor switches automatically to 1x1 layout. When the evaluation is finished the original monitor layout is restored.

Measurements are not possible for images imported from other modalities and images from external video sources (CAP images).

Selecting images for evaluation

Graphical evaluations can only be performed on a single image at a time.

- ◆ Display the required image (single selection).

Displaying the evaluation functions


The evaluation functions are available on every control panel in the **Measure Annotate** menu. If you are working on the monitor trolley, you can switch to enhanced mouse mode and open the evaluation functions from the screen menu (→ Page 171 *Evaluation in enhanced mouse mode*).

Measure Annotate menu

- ◆ Press this button.

The **Measure Annotate** menu appears.

General work steps

During graphical evaluation you use graphic objects (such as lines and text fields) that you draw in the image on the left monitor. The graphic objects are saved in the image automatically.

The following presents the basic steps you will use with graphic objects regardless of the evaluation function.

Moving the mouse pointer

Use the mouse pointer to set the position in the image where you want processing of graphic objects to take place (such as insert, select).

- ◆ Use the arrow buttons to move the mouse pointer in the image.

Vertical arrow buttons: Move up/down

Horizontal arrow buttons: Move left/right

Diagonal arrow buttons: Combined horizontal and vertical movement

– or –

Move the mouse pointer with the mouse.

Selecting

You must select an existing graphic object to process it.

- ◆ Move the mouse pointer with the arrow buttons or the mouse to the graphic object.

Highlighting of the reference points (vertices, center point) indicates that the graphic object is selected.

Panning

1 Select the graphic object.

2 Press this button.

The graphic object is linked to the mouse pointer.

3 Move the graphic object to the required position with the arrow buttons or mouse.

4 Press this button.

The link between the graphic object and the mouse pointer is removed.

Deleting (individual)

1 Select the graphic object.

2 Press this button.

The selected graphic object is deleted.

Deleting (all)

◆ Press this button.

All graphic objects in the image are deleted.

Calibration

If you perform evaluations with distance measurements, you must calibrate the image.

CAUTION

Misinterpretation of results.

Incorrect interventional decision!

- ◆ During length calibration or when measuring or interpreting distances in images, note that the acquired images represent two-dimensional projections of real objects.

Depending on the system, the calibration error is ± 0.2 mm. Under real conditions, an error of at least ± 1 mm must be assumed even if you proceed carefully (with a sufficiently large and visible calibration object, pixel-accurate marking of the end points and the calibration object and measured object in the same plane). The error is largely determined by the position of the measured object in relation to the image intensifier/detector.

Calibration object

For increased accuracy it is recommended to use a calibration object that is as large as possible, but still can be displayed completely. The edges of the calibration object must be clearly recognizable and its dimensions must be known.

CAUTION

Risk of entering the incorrect calibration. Incorrect calibration display.

Incorrect measurement of graphic objects. Incorrect interventional decision!

- ◆ The user must carefully enter the known length of the calibration object. The calibration object must be perpendicular to the X-ray beam.

Please note that the calibrating accuracy decreases if the calibration object is not located in the same plane as the measured object.

Performing calibration


- 1 Draw a distance line along the calibration object.

Be careful to mark the starting and end points with as close to pixel accuracy as possible in order to achieve optimal measuring accuracy (the line to be calibrated has to be selected/marked).

- 2 Then press this button.

The control panel switches to keyboard mode.

- 3 Enter the distance (in mm) and confirm with Enter.

A length scale is displayed in the image and is applied to the remaining images of the series.

Distances already drawn in are updated.

Measuring the distance

Using a distance line, you can measure the distance between two points in an image. Proceed carefully when setting the starting and end points in order to be as accurate as possible.

Distance drawing mode

- ◆ Press this button to activate the **Distance** drawing mode.

Setting the starting point

- ◆ Use the arrow buttons to move the mouse pointer to the starting point of the line and press this button.

– or –

Click the starting point with the left mouse button.

The starting point of the line is set.

The button changes to **Set last point**.

Setting the last point

- ◆ Use the arrow buttons to move the mouse pointer to the last point of the line and press this button.

– or –

Click the last point with the left mouse button.

The distance line is set.

The length (in mm) and the sequential number of the graphic object are displayed at the end point of the line.

The inaccuracy of a length measurement is at least ± 0.2 mm even under ideal conditions (with an optimally performed calibration and the measured object lying directly on the detector). Inaccuracy $> \pm 1$ mm must be assumed under real conditions.

5 Postprocessing

Deleting points

If the starting point or end point has been placed incorrectly, it can be deleted immediately after being drawn.

- ◆ Press this button.

The last point set is deleted from the image.

Pressing the button again deletes the starting point if available.

You can continue drawing the line.

Measuring an angle


You can define an angle by two lines, the legs of the angle that you draw on the image. Proceed carefully when setting the starting and end points in order to obtain as high an accuracy as possible.

Angle drawing mode

- ◆ Press this button to activate the **Angle** drawing mode.

Drawing an angle

- 1 Draw the first line in the image by setting the starting point and end point (→ Page 169 *Measuring the distance*).

The first leg of the angle is set.

- 2 Draw the next line in the image by setting the starting point and end point (→ Page 169 *Measuring the distance*).

The angle is set.

The drawn angle (in °) and the sequential number of the graphic object are displayed at the end points of the legs.

Either an inner or outer angle is indicated depending on the starting point and end point of the second leg.

The accuracy of the angle measurement largely depends on how carefully the user proceeds. In normal cases, the inaccuracy of the angle measurement must be assumed to be $\pm 1^\circ$. The accuracy can be improved by increasing the distances between the starting and end points of the legs drawn in.

Annotating images

You can annotate interesting or anomalous areas in the image accordingly.

Text drawing mode


- 1 Press this button to activate the **Text** drawing mode.
- 2 Move the mouse pointer to the required position for the annotation with the arrow buttons or mouse.

Entering text


- 1 Press this button.
The control panel switches to keyboard mode.

- 2 Enter the text and confirm with Enter.
The text is inserted at the position of the mouse pointer.

Inserting predefined text

Frequently used texts can simply be selected from a template once they have been created as such (→ Page 174 *Annotating images*).


- ◆ Select the required text from the selection list.
The text is inserted at the position of the mouse pointer.

Changing text

- 1 Select the graphic object with the text.
- 2 Press this button to change the text using the keyboard.
– or –

Select a new text from the selection list.

Setting the side identification

By identifying the side of the patient you avoid a mistake in image orientation. The side identification is applied automatically to the remaining images of the series.

- ◆ Switch on the **Text** drawing mode, if required.

Adding side identification

- 1 Move the mouse pointer to the required position for the identification (such as the edge of the image) with the arrow buttons or mouse.

- 2 Press the button for the correct patient side.
The selected identification **L** or **R** is inserted at the position of the mouse pointer.

Changing the side identification

- ◆ Press the button for the correct patient side.
The original side identification is deleted.
The selected identification **L** or **R** is inserted at the position of the mouse pointer.

5.2.6 Evaluation in enhanced mouse mode

When you evaluate images on the monitor trolley you can perform all evaluation functions directly on the left screen with the mouse. You only use the control panel to enter text.

5 Postprocessing

Activating enhanced mouse mode

- 1 Click this button on the left monitor (lower left).

The buttons for evaluation with the mouse are displayed (upper left).

- 2 Click the button for the required function.

Enhanced mouse mode is activated.

The control panel switches to keyboard mode.

Turning off enhanced mouse mode

- ◆ Click the button of the activated function.

– or –

End keyboard mode.

General work steps

During graphical evaluation you use graphic objects (such as lines and text fields) that you draw in the image on the left monitor. The graphic objects are saved in the image automatically.

The following presents the basic steps you will use with graphic objects regardless of the evaluation function.

Selecting You must select an existing graphic object to process it.

- 1 Click this button.

The "Selection" function is activated.

- 2 Click the graphic object.

Highlighting of the reference points (vertices, center point) indicates that the graphic object is selected.

Panning 1 Select the graphic object.

- 2 Press the mouse button and drag & drop the center point to the required position.

Changing dimensions

- 1 Select the graphic object.

- 2 Press the mouse button and drag & drop the vertices in the required direction.

Deleting (all)

- ◆ Click this button on the left monitor.

All graphic objects in the image are deleted.

Deleting (individual) 1 Select the graphic object.

- 2 Press this button on the control panel.
The selected graphic object is deleted.

Calibration


- 1 Draw a distance line along the calibration object and select it
(→ Page 173 *Measuring the distance*).

Be careful to mark the starting and end points with as close to pixel accuracy as possible in order to achieve optimal measuring accuracy.

- 2 Then click this button.

The **Calibrate** dialog box is displayed.

- 3 Enter the distance (in mm).
4 Click **OK**.

A length scale is displayed in the image and is applied to the remaining images of the series.

Distances already drawn in are updated.

Measuring the distance

Using a distance line, you can measure the distance between two points in an image. Proceed carefully when setting the starting and end points in order to obtain as high an accuracy as possible.

- 1 Click this button.
2 Click the starting point of the line.
3 Click the desired end point.

The distance line is set.

The length (in mm) and the sequential number of the graphic object are displayed at the end point of the line.

The inaccuracy of a length measurement is at least ± 0.2 mm even under ideal conditions (with an optimally performed calibration and the measured object lying directly on the detector). An inaccuracy $> \pm 1$ mm must be assumed under real conditions.

Measuring an angle

You can define an angle by two lines, the legs of the angle that you draw on the image. Proceed carefully when setting the starting and end points in order to obtain as high an accuracy as possible.

5 Postprocessing

- 1 Click this button.
- 2 Draw the first line in the image by setting the starting point and end point.
The first leg of the angle is set.
- 3 Draw the next line in the image by setting the starting point and end point.
The angle is set.
The drawn angle (in °) and the sequential number of the graphic object are displayed at the end points of the legs.
Either an inner or outer angle is indicated depending on the starting point and end point of the second leg.

The accuracy of the angle measurement largely depends on how carefully the user proceeds. In normal cases, the inaccuracy of the angle measurement must be assumed to be $\pm 1^\circ$. The accuracy can be improved by increasing the distances between the starting and end points of the legs drawn in.

Annotating images

In addition to side identification and annotating interesting or anomalous areas in an image, you can add image comments. Comments always appear centered at the bottom edge of the image and are transferred automatically to all images of the series.

- ◆ Click this button.
The input dialog box for text and comments is displayed.

Text	
L	R
xyz	
Add to predefs.	Del from predefs.
Comment	
zyx	
Add to predefs.	Del from predefs.

Setting the side identification

By identifying the side of the patient you avoid a mistake in image orientation. The side identification is applied automatically to the remaining images of the series.

- 1 Select the appropriate patient side in the text and comments dialog box.
- 2 Using the mouse, click the appropriate location in the image (e.g. at the edge of the image).
The selected identification L or R is inserted at the position of the mouse pointer.

Inserting annotations

- 1 Select the annotation from the **Text** selection list.
– or –
Click the text field of the **Text** selection list and overwrite the displayed text with the keyboard.

- 2 Click the required position in the image.

Adding a comment

- ◆ Select the text from the **Comment** selection list.

– or –

Click the text field of the **Comment** selection list and overwrite the displayed text with the keyboard.

To delete a comment from the images, use the **Backspace** button to delete it from the text field of the **Comment** selection list.

Creating predefined texts

You can save frequently used texts and comments as templates and select them from the list as needed.

- 1 Enter the text in the text field of the corresponding selection list.

Add to predefs.

- 2 Click this button.

The text is added to the list as a new entry.

Up to 20 entries each can be created in the **Text** and **Comment** selection lists.

Deleting predefined texts

- 1 Select the entry you no longer require from the corresponding selection list.

Del from predefs.

- 2 Click this button.

The text is deleted from the list.

5.2.7 Saving images

After completing postprocessing of a prior examination, all changes in the image are saved when you close the patient. If you want to save intermediate states, manually save the images in the required processing phase. The images are saved as copies of the original in the corresponding study.

To save images while an examination is in progress see (→ Page 126 *Saving and displaying images*).

Saving a single image

- ✓ Multiple images are not selected (→ Page 157 *Multiselection*).

- 1 Display the required image (single selection).

- 2 Press this button.

Saving multiple images

With the appropriate image selection you can save a portion or all images of one or more series, as well as segments of series.

Images/series

- 1 Select the required images/series in multiselection (→ Page 157 *Multiselection*).

5 Postprocessing


2 Press this button.

Partial scene

1 Show the scene.

2 Open the **Cine replay** menu.

3 Limit the replay range of the scene to the relevant section.

4 Press this button.

The images are saved from the starting point to the end point as a new scene.

5.3 Subtraction processing

When subtraction is performed, two X-ray exposures are acquired of the region of interest, one without contrast agent and one with contrast agent. Afterwards both exposures are subtracted.

In the subtraction image, areas with the same attenuation and areas with a (slightly) different attenuation, e.g. vessels with contrast agent, can be discerned clearly.

5.3.1 Postprocessing functions

The following functions are available to postprocess subtraction images during or after an examination:

- Anatomical background
- Pixel shift
- Stenosis quantification

5.3.2 Loading and displaying images

Native images from examinations acquired in Subtraction or Roadmap operating mode (optional) are stored in the local database together with the subtraction images. Both image types can be displayed simultaneously on both monitors.

Loading images

With an examination in progress, the images required for subtraction processing are already loaded and saved. For subsequent processing, load the subtraction data from the local database.

1 Select the subtraction study of the required patient in the **Patient list**.

2 Press this button in the symbol button bar on the control panel.

The selected images are loaded and displayed on the monitors.

The **PROCESSING** task card is displayed on the control panel.

Changing the monitor layout

After loading the images, 2-channel display is activated with both monitors in 1x1 layout. The left screen shows the subtraction images, the right screen the corresponding native images. This setting can be switched if necessary.

- 1 Press this button.

The **Sub / Native** menu appears.

- 2 Select the required monitor layout.

2-channel display (only for full image display on the left monitor). Simultaneous display of the subtraction image (left monitor) and corresponding native image (right monitor).

1-channel display – native images. Only native images are shown on both monitors. Applies particularly when the image overview is activated on the left monitor.

1-channel display – subtraction images. Only subtraction images are shown on both monitors. Applies particularly when the image overview is activated on the left monitor.

Displaying images

After switching the monitor to 1-channel display, use the image overview display to show the required images.

- 1 Switch on image overview if necessary.
- 2 Use the functions for image display in the same manner as for standard examinations (without subtraction).

- Changing the display mode (→ Page 154 *Changing the display mode*)

- Scrolling (→ Page 155 *Scrolling*)

In stack display the maximum fill image (Peak OP image) of the respective series is displayed.

- Reviewing a scene (→ Page 156 *Reviewing a scene*)

In 2-channel display, the subtraction scene plays on the left monitor, the corresponding native scene simultaneously on the right monitor.

5.3.3 Changing the image display

As with standard examinations (without subtraction) you can window, zoom/pan, rotate, and flip subtraction and native images (→ Page 160 *Changing the image display*).

The following presents special functions for use exclusively with subtraction data sets.

Windowing images

The change in windowing settings always refers to an image type. In 2-channel display you select whether you want to window subtraction images or native images.

- 1 Open the **Brightness Contrast** menu.

5 Postprocessing

- 2 Press the button for the required image type.
- 3 Change the window value settings using the functions in the **Brightness Contrast** menu.

– or –

- Press the button for automatic window value assignment.
If available, window values are adjusted based on the maximum fill image.

Including anatomical background (Landmark)

Normally the anatomical surroundings of vessels of interest are not visible in images that are displayed subtracted. By adding the mask image, the surrounding tissue can be highlighted more or less.

- 1 Press this button.
The **Pixel shift Landmark** menu appears.

- 2 Set the degree of inclusion with the slide control.
0%: Fully subtracted (without anatomical background)
100%: Native display

Making the image and mask coincide exactly (Pixel shift)

During acquisition of the subtraction series, the patient or the C-arm system may have moved. The images that are to be subtracted may not be congruent. The anatomical background, especially in the area of image contrasts, is therefore not canceled out.

Through a suitable shift of the mask by a few pixels, you eliminate the displacement of fill image and mask. Applying a pixel shift will also have an effect on the representation of the fill image.

The image shift with pixel shift is used in forward direction for subtraction scenes; that is, for subsequent images, not backward for previous images.
You can use manual and auto pixel shift multiple times one after the other or in combination as soon as the previously changed pixel shift is accepted.


Pixel shift Landmark menu

- 1 Display the required subtraction series (single selection).
- 2 Press this button.
The **Pixel shift Landmark** menu appears.
The subtraction series is shown in image stripe display on the left monitor.
The currently selected image of the subtraction series is shown on the right monitor.

Selecting the image

- ◆ Use the slider to display the first image with poor coverage on the right monitor.

**Automatic pixel shift**

Mark an area of obvious displacement in the subtraction image. The pixel shift needed for alignment is calculated automatically.

- 1 Press this button to switch from manual pixel shift to auto pixel shift, if necessary.

In the preview image on the right monitor, a rectangle indicates the mouse position in the center of the image.

- 2 Move the rectangle to the position in the image where misalignment is most evident.


- 3 Press this button to apply the pixel shift to the current and subsequent images.

Manual pixel shift

Shift the mask manually until the best impression of the subtraction image is achieved.

- 1 Press this button to switch from auto pixel shift to manual pixel shift, if necessary.

- 2 Move the mask with the arrow buttons until you get optimum coverage.

- 3 Press this button to apply the pixel shift to the current and subsequent images.

Peak-OP preview

Peak-OP preview allows checking the modified pixel shift before the new value is saved using **Accept pixel shift**.

- ◆ Press this button to calculate a Peak-OP image based on the current pixel shift settings.

The new Peak-OP image is displayed on the right monitor.

Undoing pixel shift

Accept pixel shift applies the changed pixel shift to the affected images and saves it. In the event of an error, you can cancel the last step or all steps.

- ◆ Press this button to cancel the last change in the pixel shift.

– or –

Press this button to restore the state prior to opening the **Pixel shift Landmark** menu.

5.3.4 Evaluating subtraction images

The following presents special evaluation functions for subtraction data sets.

Stenosis quantification

In stenosis quantification, the reduction of vascular diameter is determined through a stenosis identified in the image.

Stenosis measurements should preferably be done using subtracted images. The edge position has to be checked in native display or with anatomical background. If necessary, the vascular edges have to be corrected manually.

Opening stenosis quantification


1 Display the desired fill image (single selection).

2 Click this button on the left monitor (lower left).

The control panel switches to keyboard mode.

The **Stenosis Quantification** dialog box is displayed.

Defining the vascular course

You identify the vascular course by drawing the center line of the vessel as a polygon. With each change of direction, you set a vertex in the image. The center line of the vessel should include the section with the pathological (stenotic area) and normal vascular diameter (reference area).

1 You set the starting point and vertices with a single click of the mouse.


A line is drawn connecting the points.

2 Set the end point with a double click of the mouse.

A vascular contour is determined and drawn in the image.

Stenosis calculation

After drawing the vascular contour, the "stenosis diameter" (white line) and "reference diameter" (gray line) are determined automatically, and each is displayed as a bar cursor.

The degree of stenosis relative to the reference area is output in %:

- %D: (D = diameter) deviation of the diameter on line "Stenosis diameter" as a percentage relative to the diameter on the line "Reference diameter" $(1 - D_{sten}/D_{ref}) * 100\%$
- %A: (A = area) deviation of the calculated area of the circle as a percentage $(\frac{1}{4} * \pi * D^2)$ on the "Stenosis diameter" line, relative to the calculated area of the circle on the "Reference diameter" line $(1 - A_{sten}/A_{ref}) * 100\%$
- %G: (G = gray level) deviation of the gray level on the "Stenosis diameter" line as a percentage relative to the gray level on the "Reference diameter" line $(1 - G_{sten}/G_{ref}) * 100\%$

Please note:

- Values are limited to $\pm 100\%$
- Negative values are an indication that the lines for stenosis diameter and reference diameter are incorrectly assigned.

In addition to the imprecision caused by the parallax effect, limits to measurement precision are set that trace back to internal rounding errors. This additional imprecision generally does not exceed 10% for area measurements, 5% for diameter measurements, and 1% when measuring grayscale proportions.


With non-circular stenoses, the stenosis may only be recognizable by the gray level ratio:

	Example: Vessel cross-section (schematic)	Possible display in X-ray image (schematic)	Result of stenosis quantification
Vessel circular (no stenosis)			No stenosis found for diameter, area, gray level
Non-circular stenosis			Stenosis only found for gray level

Moving the stenosis/reference cursor


The positions of the stenosis and reference cursors can be changed individually. The stenosis is then recalculated automatically.

5 Postprocessing

- 1 Drag the slider to the desired position or click the arrow buttons.

– or –

Press these buttons on the control panel to move the stenosis cursor.

- 2 Press these buttons on the control panel to move the reference cursor.

Correcting contour lines

You change the course of the contour line by manually redrawing the vascular edge in the section to be corrected. The stenosis is then recalculated automatically after the correction.

- 1 Mark the starting point of the vascular contour to be corrected with a click of the mouse.

If the mouse pointer is located over a valid vascular position, it takes the shape of a cross.

- 2 Mark the end point of the vascular contour to be corrected by double clicking. The vascular contour between these points is changed.

Moving the stenosis text block

After the stenosis is drawn in the image, you can move the text block to another position, for example, if it is covering an important detail.

- ◆ Click the text block; drag the text block to the new position, and click again at the new position.

A number indicates the stenosis cursor to which the text block belongs.

Confirming stenosis quantification

- ◆ Click **OK** in the **Stenosis Quantification** dialog box.

If necessary you can evaluate additional stenoses in the image in the same manner.

Deleting stenosis calculation (individual)

- 1 Click the stenosis to be deleted. The stenosis is marked by small squares.
- 2 Click this button in the **Stenosis Quantification** dialog box. The selected stenosis graphic with evaluation is deleted. Other stenosis graphics are retained.

Saving the stenosis image

- ◆ Click this button in the **Stenosis Quantification** dialog box. The stenosis graphics are “burned into” the image. The image is stored as a new image and is assigned to the study.

The results of the stenosis quantification tool highly depends on the image content. So the results are influenced by the body anatomy (thickness and vessel form and routing), the projection (angle, distance to focus, etc.), the achieved image contrast (dose), the contrast medium and many other factors.

The displayed values from the stenosis tool always need to be checked for plausibility. This can be achieved by Subtraction images from a different projection angle and additional indicators like (EEG, ECC).

6 Documentation

6.1	Filming/Printing	186
6.1.1	Selecting images	186
	Single image	186
	Multiselection	186
6.1.2	Printing using default settings	187
	Starting printing	187
6.1.3	Printing with changed settings	187
	Input dialog box for print settings	188
	Changing print settings	188
	Displaying a print preview	189
	Starting printing	189
6.2	Exporting	189
6.2.1	Export procedures	190
	Standard export	190
	Changing the export settings	190
	Automatic export	190
6.2.2	Export destinations	190
	Removable device	190
	Other workstation	191
	Archive	191
6.2.3	Overview of transfer capabilities	191
6.2.4	Selecting images	192
	Images from an examination in progress or postprocessing	193
	Complete studies	193
6.2.5	Preparing the removable devices	194
	Working with data media	194
	Drives	194
	Inserting a CD/DVD	194
	Connecting the USB storage medium	194
6.2.6	Exporting using default settings	194
	Starting export	194
6.2.7	Exporting with changed settings	195
	"Export to" dialog box	195
	Changing the export settings	196
	Starting export	198
	Establishing default settings	198
6.3	Checking the data transfer	199
6.3.1	Indicators in the status bar	199
6.3.2	Viewing and processing jobs	199
	Job status list	199
	Influencing job performance	200
6.4	Reports	201
6.4.1	Selecting a study	201
6.4.2	Radiation Summary Report	202
	Opening the report	202
	Exporting the report	202
	Printing a report	203

6.1 Filming/Printing

You can expose the images of a study on film or print them on paper for documentation and reporting.

Make sure the printer is switched on before you send images to print.

Filming and printing are the same process except that they use different output devices. If the term "filming" or "printing" is used alone in the following, the description applies equally for the other process.

Standard printing

You can print the displayed image or a selection of multiple images simply by pressing a button. The print job is processed with standard settings (simplified print procedure).

Changing print settings

If you want to check or change settings prior to printing, use the "enhanced print procedure". You can display a print preview and adapt the print layout, for example. In a second step, you forward the print job to the required printer.

6.1.1 Selecting images

Before you start a print job, you have to select the images.

You have access to the currently loaded images through the **PROCESSING** task card on the control panel. Single and multiple image selections are possible (→ Page 157 *Selecting images*).

- ◆ Select the suitable display mode in the **Overview** menu:

Image stripe display if you want to select multiple images in the current series.

Stack display if you want to select one or more series.

Single image

- ◆ Select the required image or series in the **Overview** menu.

If the left monitor has been changed to full image mode, the image displayed there is selected.

Multiselection

- ◆ Press this button in each case to add the image or series selected in the image overview to the multiple selection.

– or –

Switch on multiselection in the **Overview** menu.

Select the images or series.

Images of a multiselection are marked on both monitors by a white dot.

The **Multiple Selection** menu has to be open/active during printing; otherwise, only the currently marked image/scene is printed.

6.1.2 Printing using default settings

The following default settings are used with the simplified print procedure:

- Page size and layout (page partitioning) according to the configuration settings for the **Default Printer** (→ Page 213 *DICOM Properties - Printing*).
- Handling of subtraction images and image text display as in the last print job. These settings can be changed in the enhanced print procedure (→ Page 188 *Changing print settings*).

Starting printing

- ◆ Press this button.

The print job is started.

In the **Patient list** the selected images are marked as "printing" ("p" flag).

When the print job is complete the images are marked as "printed" ("P" flag).

After printing, the images are released for automatic deletion if necessary. Depending on the hard disc fill level, the images are deleted on the specified date (→ Page 218 *Deleting images automatically*).

The "printed" flag is set as soon as the images are successfully transferred to the printer driver. Not all printers (e.g., paper printers) can solve printing problems themselves. The image printout may be lost! Verify that the printouts are available before you delete images.

6.1.3 Printing with changed settings

In the enhanced print procedure you can change the following settings:

- Printer selection
- Layout settings
- Image text display
- Handling subtraction images
- Number of copies

Input dialog box for print settings

(1) Print settings

(2) Buttons for print, print preview, and print status

◆ Click this button on the left monitor (lower left).

The input dialog box for print settings is displayed on the monitor.

Changing print settings

The input dialog box for print settings initially contains the settings used for the last print session.

Changes to the settings for subtraction images and image text display are used as new default settings for standard printing.

Printer The **Printer** selection list contains the paper and DICOM printers configured in the system.

◆ Select the required printer.

Layout settings Individual print layouts under an individual name are created in the configuration for each printer (→ Page 216 *Print Layout task card*).

◆ Select the layout settings in the selection list.

Name: Stored print layouts for the selected printer

Size: Available page sizes for the selected printer

Layout: Available page partitions for the selected printer

Image text display The scope of image text information is defined in the configuration for each print layout (→ Page 216 *Image text settings*).

Anonymous

1 Select the **Anonymous** option if you want to avoid assigning printed images to a patient.


"Anonymous" appears in the hardcopy instead of the patient name. The remaining demographic data are hidden.

2 Select the **Hide Text** option if all image text is to be hidden in the hardcopy.


Subtraction images

Subtraction images are usually printed in native display.

◆ If needed, select the subtracted display for printing.

Number of copies

◆ Define the number of copies.

Up to 9 copies per print job are possible.

Displaying a print preview

In the print preview you can see how the images will look later on the exposed film or paper hardcopy. You can also change the position of individual images where necessary.

1 Click this button in the input dialog box for print settings.

The images appear on the left monitor in the selected print layout.

2 Use the mouse wheel to display the individual pages of the print job.

3 Use the mouse's drag and drop function to arrange the images in the selected page layout.

Starting printing

◆ Click this button in the input dialog box for print settings.

The print job is started.

In the **Patient list** the selected images are marked as "printing" ("p" flag).

When the print job is complete the images are marked as "printed" ("P" flag).

After printing, the images are released for automatic deletion if necessary. Depending on the hard disk fill level, the images are deleted on the specified date (→ Page 218 *Deleting images automatically*).

The "printed" flag is set as soon as the images are successfully transferred to the printer driver. Not all printers (e.g., paper printers) can solve printing problems themselves. The image printout may be lost! Verify that the printouts are available before you delete images.

Print status

You can track and influence job processing for print jobs on DICOM printers.

◆ Click this button in the input dialog box for print settings.

The **Job status list** opens (→ Page 199 *Job status list*).

6.2 Exporting

After an examination or postprocessing, the saved images are stored on the hard drive (local database).

This section explains how to save images and patient data from the local database, send them within the network, and save them to data media.

6.2.1 Export procedures

Standard export

During an examination or postprocessing, you can export the displayed image or a selection of multiple images simply by pressing a button. The export job is processed using default settings (simplified export procedure).

Changing the export settings

If you want to check or change settings prior to export, use the "enhanced export procedure". You can then change the export destination and data format.

Automatic export

Depending on the configuration setting, at the end of the examination acquired images are sent automatically to a default address (→ Page 213 *DICOM Properties - Sending*).

The address and export settings used are the same as for default export.

6.2.2 Export destinations

Depending on the job setting, you select one of the following export destinations in the enhanced export procedure.

Removable device

You write data to removable devices to save it short term or pass it on.

The following media types are supported for saving data:

- DVD-R (minus R), DVD+R
- CD
- USB storage media (external hard drives, sticks)

The accidental use of other media, e.g., DVD-RW, can result in an error message or impair functionality.

CAUTION

Risk that CD/DVD-R media used for exporting can become damaged or be unreadable on other CD/DVD drives.

Data loss or apparent data loss.

- ◆ Use removable devices as transfer media only. Only data media approved by Siemens may be used.

In general, do not use storage media that require their own power supply.
In addition, do not use USB storage media containing an autorun.inf file.

As a rule, you cannot release radiation during the CD/DVD write process. In exceptional cases (emergencies), you can release radiation in fluoroscopy operating mode. However, radiation may be disrupted.

Moreover, the export process of subtracted images to CD/DVD is interrupted during radiation release and other processing-intensive procedures (e.g., play LSH). In those cases images are stored as raw data.

You are therefore urgently advised to start write processes outside of examination hours (e.g. at the end of the office day).

Other workstation

If your system is connected to a network, you can send patient and examination data to other workstations via the network. The data can be processed or used for diagnosis there.

Archive

To archive, you send the selected patient and examination data to a DICOM network node that has been established as an archive. You can import archived data back to your local database whenever you need them.

CAUTION

Use of the imaging system hard drive as a long-term archive.

Loss of image data.

- ◆ Do not use the imaging system hard drive as a long-term archive. Follow the local regulations regarding the archiving of X-ray images, while taking into account the use of DICOM archive nodes.

CAUTION


Loss of image data.

Risk of unnecessary radiation exposure!

- ◆ The hard disk of the imaging system is not suitable for long-term archiving of patient and image data. There are statutory requirements governing the archiving period, data availability, and data security (data integrity), as well as recommendations concerning fire protection or water damage for the archiving of image data. The operator of the archive is responsible for observing these regulations.
- ◆ Because of advances in technical development and statutory requirements, the storage of image data and access to it cannot be realized with a single storage and media technology alone. Therefore, the migration of data may be necessary under the responsibility of the operator of the digital archive.

6.2.3 Overview of transfer capabilities

The figure below shows the data backup and transfer options available to you.

- (1) Export to DICOM network nodes (Archive or other workstation)
- (2) Export to removable devices (CD/DVD, USB drive)
- (3) Import from DICOM network nodes (Query & Retrieve) (→ Page 150 *Data in the network*)
- (4) Importing from removable devices (→ Page 149 *Data from removable devices*)

Please remember that not all transfer options may be available on your system. The devices and network nodes available depend on the individual configuration of your system and the options installed.

6.2.4 Selecting images

Before you start an export job, you have to select the images.

You can only export objects that are stored in the local database. If data are to be transferred from one data medium to another, they must first be imported.

Exporting data from other modalities is not possible however.

Images from an examination in progress or postprocessing

You have access to the currently loaded images through the **PROCESSING** task card on the control panel. Single and multiple image selections are possible (→ Page 157 *Selecting images*).

- ◆ Select the suitable display mode in the **Overview** menu:

Image stripe display if you want to select multiple images in the current series.

Stack display if you want to select one or more series.

Single image

- ◆ Select the required image or series in the **Overview** menu.

If the left monitor has been changed to full image mode, the image displayed there is selected.

Multiselection

- 1 Switch on multiselection in the **Overview** menu.

- 2 Select the images or series.

Images of a multiselection are marked on both monitors by a white dot.

– or –

Press this button to add the image or series selected in the image overview to the multiple selection.

The **Multiple Selection** menu has to be open/active during export; otherwise, only the currently marked image/scene is exported.

Complete studies

To save or archive complete studies, select the data sets in the **Patient list** (→ Page 143 *Searching for and selecting patient data*).

- 1 Change to the **PREPARATION** task card if necessary.

- 2 Open the **Patient list**.
- 3 Search for the patient, sorting and filtering the list as necessary.
- 4 Select the required study or studies.

6.2.5 Preparing the removable devices

Depending on the routine you use, one or several drives and connections for external data media will be established for your system.

Working with data media

Follow the instructions of the manufacturer for handling and storage of CD/DVDs.

Drives

Your system comes with a CD/DVD writer for data storage.

Inserting a CD/DVD

CDs/DVDs can be written to on one side only.

- ◆ Place the CD/DVD into the drive with the label facing upwards.

Never shut down the Cios Fusion or disconnect the monitor trolley from the C-arm system while data is being written to CD/DVD.

Connecting the USB storage medium

CAUTION

Risk of connecting/disconnecting USB devices during the exposure.

Such behavior can lead to system instability. This can also impact other processes, such as exposure jobs in progress. Loss of data, interference in exposure process.

- ◆ Do not insert or disconnect USB devices during exposure or other critical processes.

CAUTION

Risk of disconnecting USB devices during ongoing data transfer.

Loss of data and damage to the data medium possible.

- ◆ Do not disconnect the USB device during active data transfer (such as export).
- ◆ Use removable devices as transfer media only. Removable devices are not permitted for archiving purposes.

- ◆ Connect the USB plug into the USB jack.

6.2.6 Exporting using default settings

During the simplified export procedure, the settings established as defaults are used. The default settings can be changed in the enhanced export procedure (→ Page 196 *Changing the export settings*).

Starting export

On the control panel there is an export button whose appearance depends on the default export destination set.

- ◆ Press this appropriate key.

The export job is started.

After successful completion of the export job to a network node, the selected images are marked as sent ("S" flag) or archived ("A" flag) in the **Patient list**.

After exporting, the images are released for automatic deletion if necessary. Depending on the hard disc fill level, the images are deleted on the specified date (→ Page 218 *Deleting images automatically*).

The "sent" or "archived" flag is set as soon as the images are successfully transferred to the network nodes. Even when the addressee works with Storage Commitment, the flags only identify the receipt and storage of the data on the recipient's hard drive. A misinterpretation of the flags can result in a loss of data during the prescribed storage period. For this reason, follow the regulatory requirements regarding the archiving procedure.

Note that the exported images are not transmitted as raw data, but including all of the applied processing steps. They cannot be restored to their original state. Do not use these images for the primary diagnosis!

6.2.7 Exporting with changed settings

In the enhanced export procedure you can change the following settings:


- Export destination selection
- Data format
- Image properties
- Handling subtraction images
- Saving with a DICOM Viewer

"Export to" dialog box

- ◆ Click this button on the left monitor (lower left).

The **Export to...** dialog box appears on the monitor.

- (1) Export settings
- (2) Storage capacity or removable devices
- (3) Buttons for export, export status, and close

Changing the export settings

In the **Export to...** dialog box you can first enter settings that will be used as the defaults.

- Export destination** ♦ Select the required destination medium from the **Target Type** selection list.
- DICOM node:** Another workstation, archive, PACS
 - CD/DVD:** Inserted CD or DVD
 - USB device:** Connected USB data medium

When DVD recording is active, export to DVD is not possible.

- Export address** When selecting "DICOM node" you have to enter a network address configured in the system as the destination.

- ♦ Select the network node from the **Target** selection list.

- Data format** The systems and applications in which the exported data can be opened and processed depend on the data format.

- ♦ Select the data format in the **Export Format** selection list.
- Private:** Format that is only readable by FLC imaging systems
- DICOM:** Standardized file format in medical imaging

TIFF: PC graphics format

AVI: PC video format

When "DICOM node" is selected as the export destination, only the DICOM format can be used.

Image properties

Depending on the data format, prior to export images can be prepared differently in terms of graphics and text embedding as well as file size.

When **Private** is selected as the data format, the data is exported only in its original state.

Export images as:

- Processed ☐
- Anonymized ☐
- Downsized ☐

- ◆ Select the image properties in the **Export images as** area.

Processed: Changes in the image display compared to the original are used (such as window values, image rotation/flip, zoom, pan, edge filter, set black function) and evaluation graphics are "burned" into the image.

Anonymized: "Anonymous" appears in the image instead of the patient name. The remaining demographic data are hidden.

Downsized: The files are made smaller by reducing the image resolution 50% and the bit depth to 8 bits/pixel.

To save the image brightness and contrast, select **Processed**.

Subtraction images

You can export images of a specific image type or multiple image types from subtraction series.

Export subtracted images as:

- Subtracted ☒
- Native ☐

- ◆ Select the image types in the **Export subtracted images as** area.

Subtracted: Images in subtracted display

Native: Native images without mask and without Peak OP image

Radiation Summary Report

Two different file formats are possible if the radiation summary report is exported in DICOM format.

Export Dose Report as:

- DICOM dose report ☒
- DICOM SC image ☐

- ◆ Select the file format in the **Export Dose Report as** area.

DICOM dose report: Structured Report (DICOM SR)

DICOM SC image: Image with "burned in" text


DICOM Viewer

When exporting to a removable device, a DICOM Viewer is saved to the CD/DVD or USB medium together with the image data. This allows for viewing the exported images on any computer. The DICOM Viewer starts directly from the data medium (automatically after inserting/connecting). No installation of files takes place on the computer in question.

Provide DICOM viewer ☐

- ◆ Deselect this option if you do not want to save the DICOM Viewer.

Scope of data


- ◆ Specify the scope of the images to be exported.

Selected: all selected images (especially individually selected images)

Marked: all marked images (with white dot)

All: all images in the study

Checking disk space

If a removable device is indicated as the export destination, the current storage capacity and required disk space are displayed.

- 1 Use the display to check where there is sufficient space available on the medium.
- 2 Change the data medium if necessary before starting the export.

Starting export

- ◆ Click this button in the **Export to...** dialog box.

The export job is started.

After successful completion of the export job to a network node, the selected images are marked as sent ("S" flag) or archived ("A" flag) in the **Patient list**.

After exporting, the images are released for automatic deletion if necessary. Depending on the hard disk fill level, the images are deleted on the specified date (→ Page 218 *Deleting images automatically*).

The "sent" or "archived" flag is set as soon as the images are successfully transferred to the network nodes. Even when the addressee works with Storage Commitment, the flags only identify the receipt and storage of the data on the recipient's hard drive. A misinterpretation of the flags can result in a loss of data during the prescribed storage period. For this reason, follow the regulatory requirements regarding the archiving procedure.

Export status

You can track and influence job processing for all export jobs.


- ◆ Click this button in the **Export to...** dialog box.

The **Job status list** opens (→ Page 199 *Job status list*).

Establishing default settings

Export jobs using the simplified procedure (standard export), automatic export, and sending the radiation summary report at the end of the examination, are all performed using default settings.

- 1 Select the required export settings in the **Export to...** dialog box.
- 2 Select the **Set as Default** option.

- 3 Click this button.

Changed default settings are stored.

6.3 Checking the data transfer

All jobs for importing, exporting, and printing/filming data are performed in the background.

The processing of jobs can be viewed in the status bar, enabling you to react in case of error.

6.3.1 Indicators in the status bar

The corresponding icon is displayed in the status bar on the left monitor (lower left) for data transfers in progress. The icons indicate the action being performed and whether an error has occurred.

Action	Not active	Active	Error
Write to CD/DVD	No icon		No icon
Read from CD/DVD	No icon		No icon
Write to USB device	No icon		No icon
Read from USB device	No icon		No icon
Export to DICOM network node (sending)	No icon		
Import from DICOM network nodes (receiving)	No icon		No icon
Printing/filming to DICOM printer	No icon		

6.3.2 Viewing and processing jobs

In the **Job status list** you can learn about the following jobs running in the background – sorted by transfer type:

- Printing/filming to DICOM printer (DICOM Film)
- Export to DICOM network node (DICOM Send)
- Import from DICOM network nodes (DICOM Retrieve)
- Write to CD/DVD
- Read from CD/DVD
- Write to USB storage medium
- Read from USB storage medium

Job status list


Several methods are available for opening the **Job status list**:

- ◆ Click the **Status** button in the **Import from**, **Export to** dialog box, or the print settings dialog.

– or –

Click the icon for the applicable data transfer on the left monitor.

The **Job status list** is displayed. The task card for the applicable data transfer is displayed.

- (1) Tabs for selecting the transfer type
- (2) Details and status of the transfer jobs
- (3) Area to show the course of job processing
- (4) Buttons for deleting the list or selection and to close

Influencing job performance

The work steps available in the list depend on the transfer type and the status of a job.

Selecting jobs

Some actions are limited to the transfer jobs selected in the list.

- ◆ Click the required jobs in the list (press and hold **Ctrl** for multiple selection).

Restarting DICOM Send

Failed export jobs to DICOM network nodes can be repeated with the same settings or redirected to a new address.

- ◆ Click **Send**.

The selected jobs are restarted with the original destination address.

– or –

Click **Redirect**.

The selected jobs are diverted to the DICOM network node currently set in the **Export to...** dialog box.

Aborting jobs ♦ Click **Abort**.

All transfer jobs are cancelled and the start of other planned jobs in the list is prevented.

A data transfer in progress for which there are additional images (such as a print job with multiple images on a page) is not cancelled.

Clearing lists ♦ Click **Clear List**.

All jobs that have ended or for which data transfer preparation has been completed are hidden.

– or –

Click **Clear selected**.

Of the **selected** jobs in the list, all that have ended or for which data transfer preparation has been completed are hidden.

6.4 Reports

Certain examination data are saved in the form of structured reports. Cios Fusion offers the following types of reports:

■ Radiation Summary Report

Use: The reports serve to document examination and treatment data. They are mostly generated automatically. Depending on the type of report you can read, edit/supplement, print, and export them.

Prerequisite: The examination of a registered patient is the prerequisite for a report to exist.

Formats: Reports are saved in two different formats. This allows them to be opened in different applications.

■ Secondary Captures (DICOM SC)

Allows loading into a viewer, for example; PACS compatible.

■ Structured Report (DICOM Dose)

Data format that is compatible with other DICOM systems; reports must be exported to these systems in DICOM-SR format.

6.4.1 Selecting a study

One radiation summary report is created in the local database for each study.

♦ Select an image or series of the required study in the **PROCESSING** task card on the control panel.

– or –

In the **Patient list** select the study if the applicable data have not been loaded.

6.4.2 Radiation Summary Report

The radiation summary report contains the cumulated values for the number of exposures, duration, and dose over the course of an examination. The generation of this report is a fully automated process.

(1) Patient Name 141, Emergency
Date of Birth 17-02-2017
Patient ID E_141
Sex O
Ward
Request ID
Accession No.
Study Description Intervention Id
Perf. Physician
Institution Name
Number of series/images 25 / 7138
Accumulated Values
Fluoro time 00:10:14
Dose Area Prod. ($\mu\text{Gy}\cdot\text{m}^2$) 2166.98
Ref. Air Kerma (mGy) 51.50

(2)

KV	mAs	X-Ray length ...	Pulse rate (p/s)	Dose Le...	Scene I...	Format	Char. Curve	Date and Time	
80	0.30	5.3		0.548		Mag 0	DR-CARD ...	17-02-2017 14:16:32	K4
80	0.61	4.9		0.548		Mag 0	DR-CARD ...	17-02-2017 14:17:03	K4
80	1.23	4.7		0.548		Mag 0	DR-CARD ...	17-02-2017 14:17:16	K4
79.9	1.95	7.9		0.548		Mag 0	DR-CARD ...	17-02-2017 14:17:29	K4
78.9	0.04	5.2	30.0	0.012	2.5	Mag 0	FL-CARD HP	17-02-2017 14:17:40	K4
60.5	0.13	5.3	15.0	0.012	6.3	Mag 0	FL-HC	17-02-2017 14:23:49	St
72.2	2.11	21.1		0.548		Mag 0	DR-HC	17-02-2017 14:25:23	St
72.7	2.42	24.6		0.548		Mag 0	DR-HC	17-02-2017 14:25:46	St
72.7	2.44	24.8		0.548		Mag 0	DR-HC	17-02-2017 14:25:53	St
61.1	0.13	5.5	15.0	0.012	0.3	Mag 0	FL-HC	17-02-2017 14:28:18	St
60.5	0.12	5.0	15.0	0.012	7.7	Mag 0	FL-HC	17-02-2017 14:28:30	St
72.1	2.12	21.1		0.548		Mag 0	DR-HC	17-02-2017 14:28:39	St
72.7	2.41	24.6		0.548		Mag 0	DR-HC	17-02-2017 14:28:44	St
60.5	0.12	4.9	15.0	0.012	2.2	Mag 0	FL-HC	17-02-2017 14:31:55	St
72.1	2.11	21.1		0.548		Mag 0	DR-HC	17-02-2017 14:33:20	St
60.5	0.12	5.0	15.0	0.012	7.3	Mag 0	FL-HC	17-02-2017 14:33:26	St
60.5	0.12	4.9	15.0	0.012	2.5	Mag 0	FL-HC	17-02-2017 14:34:47	St
69.4	0.40	11.4	15.0	0.012	6.2	Mag 0	FL-HC	17-02-2017 14:36:30	St
125	0.30	11.4	15.0	0.012	2.7	Mag 0	FL-HC	17-02-2017 14:36:44	St
60.8	0.14	4.9	15.0	0.012	6.1	Mag 0	FL-HC	17-02-2017 14:36:59	St
46.6	0.06	8.3	15.0	0.012	3.2	Mag 0	FL-HC	17-02-2017 14:46:33	St
69.7	0.53	11.3	15.0	0.012	0.5	Mag 0	FL-HC	17-02-2017 14:46:54	St
69.7	0.56	11.4	15.0	0.012	7.3	Mag 0	FL-HC	17-02-2017 14:46:58	St
69.7	0.53	11.3	15.0	0.012	12.2	Mag 0	FL-HC	17-02-2017 14:47:24	St
49.8	0.07	7.4	15.0	0.012	6.8	Mag 0	FL-HC	17-02-2017 14:47:49	St
69.9	0.55	11.8	15.0	0.012	4.3	Mag 0	FL-HC	17-02-2017 14:48:09	St
69.9	0.55	11.8	15.0	0.012	4.5	Mag 0	FL-HC	17-02-2017 14:48:15	St
49.6	0.09	7.3	15.0	0.012	4.4	Mag 0	FL-HC	17-02-2017 14:48:52	St
51.3	0.03	7.4	15.0	0.012	3.8	Mag 0	FL-HC	17-02-2017 14:49:04	St

(3) ☒ Print detailed dose report
☐ Print on first sheet
Export Print Close

- (1) Patient and examination data, cumulative radiation data
- (2) Detailed information on the acquired series and images
- (3) Buttons for scrolling in the report, exporting, and printing the report, as well as closing the dialog box

Opening the report

- ✓ A study or a series/image is selected.

- ◆ Click this button on the left monitor (lower left).

The radiation summary report is displayed on the monitor.

Exporting the report

The default settings established in the **Export to...** dialog box are used when exporting the report. This is particularly true of the file format when DICOM format is set.

Export


- ◆ Click this button.

The report is exported to the default address.


Printing a report

When you send the report directly to the default printer, the settings established for standard printing are used. To use other printer settings or a different printer, open the print preview first.

Standard printing

☐ Print on filmsheet


- 1 Deselect the **Print on filmsheet** option, if necessary.

Print

- 2 Click this button.

The report is printed on the default printer.

Changing print settings

☒ Print on filmsheet

- 1 Select the **Print on filmsheet** option.

Print

- 2 Click this button.

The report is opened in print preview. Make the necessary print settings and start the print job here (→ Page 188 *Changing print settings*).

7 Configuration

7.1	User Settings	207
7.1.1	Opening/closing a configuration window	207
7.1.2	Site Info	207
	Institution Name	208
	Institutional Department Name	208
	Hospital Address	208
	Patient Demographic Data	208
	Daylight saving time	208
	Application selection	209
7.1.3	Patient registration	209
	Name-display format	209
7.1.4	DICOM Properties	210
	DICOM Properties - General	211
	DICOM Properties - Protocol codes	212
	DICOM Properties - Sending	213
	DICOM Properties - Printing	213
7.1.5	Examination properties	214
	Set Black	214
	Nominal Power Shot	214
7.1.6	Image display and print properties	214
	Display Layout task card	215
	Print Layout task card	216
	Image text settings	216
7.1.7	Deleting images automatically	218
	Auto delete	218
7.1.8	Password for radiation release	219
7.2	Examination settings (PEX Editor)	219
7.2.1	Introduction	219
	PEX Editor	220
	PEX database	220
	Online/offline	220
	Data categories	220
	Names	220
	Storage capacity	221
	Access control	221
	Access modes	221
	Import/Export	221
	Automatic backup	221
7.2.2	Starting the PEX Editor	221
	PEX Editor layout	223
	General work steps in the PEX Editor	223
7.2.3	Editing application groups	224
	Assigning applications to application groups	225
	Removing applications from an application group	225
	Establishing the default application group	225
7.2.4	Editing applications	225
	Configure triplets for applications	226
	Defining other parameters	226
7.2.5	Parameter for fluoroscopy and exposure	227
	Parameter	227
7.2.6	Defining parameter sets	228

7 Configuration

7.2.7	PEX service functions	229
	Configuration parameters	229
	Changing the PEX Editor password	229
	Backing up and restoring the PEX database	230

7.1 User Settings

7.1.1 Opening/closing a configuration window

PREPARATION

1 Change to the **PREPARATION** task card on the control panel if necessary.

2 Press this button.


The control panel switches to keyboard mode.

Manage Tool is displayed on the monitor.

3 Click this icon on the **Settings** subtask card in the **Manage Tool**.

The menu for defining user settings is displayed.

4 Double-click the required settings icon to open the respective settings windows.

5 Change the settings as required. See following pages.

6 Click **Apply** to accept the changes.

7 Click **OK** to confirm the changes and close the settings window.

8 Click **Cancel** to reject your changes and close the settings window.

9 Press this button again to exit the **Manage Tool**.

7.1.2 Site Info

Site Info

In the **Site Info** dialog box you establish site-specific settings, such as the name and address of your facility.

7 Configuration

Site Info

Institution Name: Institute

Institutional Department Name:

Hospital Address:

Patient Demographic Data:

- ☒ Accession No.
- ☐ Request ID
- ☐ Patient ID
- ☐ Date of birth

Daylight saving time:

- ☒ Automatic activation

Application selection:

- ☐ Select default application after patient close

OK Apply Cancel

Institution Name

- ◆ Enter the name of your hospital or practice (64 characters maximum).

A change of the **Institution Name** will have no effect on already acquired images.

Institutional Department Name

- ◆ Enter the name of your department or ward (64 characters maximum).

Hospital Address

- ◆ Enter the address of your hospital or practice.

Patient Demographic Data

- ◆ Select the patient data to display on the monitor in addition to the patient name:

Accession No.

Request ID

Patient ID

Date of birth

Daylight saving time

- ◆ Select **Automatic activation**.

When **Automatic activation** is selected, the system automatically switches to daylight saving time.

Application selection

- ◆ Activate the **Select default application after patient close** check box.

When **Select default application after patient close** is selected, the system automatically switches to the default application after closing patient.

7.1.3 Patient registration

In the **Patient Registration Input Fields** dialog box you establish what information can be entered during patient registration.

Entries can only be made with the mouse.

- ◆ Activate or deactivate the respective entries.

Visible: Entry will be shown in the **Data Entry Dialog** window.


Mandatory: Entry must be completed in the **Data Entry Dialog** window. (Entry will be shown in the **Data Entry Dialog** window in **bold letters**).

Name-display format


Sets the parts of the patient name to be displayed:

- 1 Click **Change**.

7 Configuration

- 2 Select the parts of the patient name that should be displayed (maximum five).
- 3 Click the corresponding check box if a comma should be placed between the parts of the name.

- 4 Click **OK**.

The patient name composition that has been established is displayed.


This composition will be used when displaying patient names in images and on hardcopies/filmsheets.

7.1.4 DICOM Properties

In the **DICOM Properties** dialog box you establish the settings for DICOM functions, such as saving, sending, and printing of images.

DICOM Properties - General

**Storage Commitment**

Waiting time: The time the imaging system waits for storage commitment after a study/images have been sent to an archive server.

Number of retries: If the archive server is not available, the imaging system performs the selected number of retries.

If a DICOM node for which storage commitment is configured does not send a confirmation within the configured waiting time/number of retries, neither the **S** (sent) or **A** (archived) flag is set.

Worklist query

RIS Timer: You can set the update interval (in minutes) for receiving worklist data from the HIS/RIS (hospital/radiology information system) by specifying the update interval (in minutes).

Possible settings: Update: 0 - 1440 minutes (24 hours), 0 = no update

Modalities to query: Selection of modality/modalities to use for worklist queries.

Any modality with a two-letter name can be entered for the worklist query (**Other**).

Detailed dose information in MPPS

Send exposure dose sequence: The exposure dose sequence is documented in the performance documentation (MPPS).

Extended dose info.: The dose area product and the name of the application are listed for every exposure.

The MPPS settings can only be changed in the Service Settings and not in this dialog box.

7 Configuration

DICOM Properties - Protocol codes

DICOM Properties

General Protocol codes Sending Printing

Edit Action Item

Protocol Code 0008,0100 Code V

Code Value 0008,0100 Code V

Code Meaning 0040,0008>0008,0104 Code Meaning

Code Designator 0008,0102 Design

Code Version 0008,0103 Versio

Remove Insert

OK Apply Cancel

For MPPS Close, you can define protocol codes.


Adding a new protocol code

- 1 Enter the code: **Code Value**, **Code Meaning**, **Code Designator**, **Code Version**.
 - 2 Then click **Insert**.
- A new **Protocol Code** is created.

Deleting a protocol code

- ◆ Select the entry to be deleted from the **Protocol Code** list and then click **Remove**.

DICOM Properties - Sending

Images for autosending


- ◆ Select which images are to be sent automatically:

Off: No images will be sent automatically.

All: All acquired images will be sent automatically.

Marked: Only marked images will be sent automatically.

DICOM Properties - Printing


- Default Printer**
- **Target:** Select one of the configured printers
 - **Layout:** Select one of the currently defined DICOM layouts
(Specification: First layout of the default printer)

7.1.5 Examination properties

Examination
Properties

Activate special performance tests to check your X-ray system in the **Examination Properties** dialog box.

Set Black

Offset: You can set an offset which will be taken into account when applying the digital collimator function.

Digital collimator function: The areas outside the collimated areas are automatically blackened if so configured for the selected application.
(→ Page 225 *Editing applications*)

Nominal Power Shot

Enable Nominal Power Shot: For VA (Veteran's Health Administration) hospitals in the United States it is required to check the nominal power of an X-ray exposure at 100 kV for 100 ms in accordance with the IEC. When this function is active, the generator creates an exposure with the nominal power when triggering a single image exposure. The function is deactivated again automatically after triggering an exposure.

To activate the function, proceed as follows:

- The check mark for **Enable Nominal Power Shot** has to be set in the user settings.
- A patient has to be registered.
- A new application has to be selected (regardless of which).
 - The Nominal Power Shot is active one time. (Single image and 100 kV are displayed at this time).

7.1.6 Image display and print properties

Image Display
and Print Pr...

In the **Image Display and Print Properties** dialog box you establish the information to be displayed in the images when shown on the monitor and in the hardcopy.

Display Layout task card


Image layout **Display image text:** Show the selected text elements in the image on the monitor.

7 Configuration

Print Layout task card

- Image layout**
- **Name:** Print layout name
 - **Type:** Type of print layout (Paper, DICOM)
 - **Size:** Size of the paper/filmsheet in inches
 - **Layout:** Page partitioning
 - **Enlarge patient name:** Print the patient name larger than the other text

Image text settings

- Top left** Patient information:
- **Patient name:** Patient's last name, first name
 - **Patient ID:** Patient's identification number
 - **Date of birth:** Patient's date of birth
 - **Ward:** Ward of examination
 - **Request ID**
 - **Accession No.**

Bottom left Exposure information:

- **Dose area product:** Dose area product
- **kV:** Exposure kV
- **mA / mAs:** Exposure mA/mAs
- **Dose value / Sensitivity class:** Dose value or sensitivity
- **II/FD image size:** Input format/zoom level
- **Image no.:** Image number
- **Application:** Application program name

Bottom center ■ **Comment:** Image comment**Top right** System information:

- **Institution name:** (→ Page 207 *Site Info*)
- **System Type:** Product name
- **Physician 1 / Physician 2:** Examining physicians
- **Date of exposure:** Date acquired
- **Time of exposure:** Examination date and time

Bottom right Image information:

- **Matrix size:** Image matrix
- **Image type:** Image type
- **Mask Number:** Number of the mask image
- **LUT level:** Name of the look-up table (LUT)
- **Edge level:** Edge enhancement module name
- **Window values:** Window values
 - **WW:** Window width
 - **WC:** Window center
- **Window values for subtracted images:**
 - **WC:** Window contrast
 - **WB:** Window brightness
- **Pixel shift values:** Pixel shift values X: .../Y: ...

Labeling

Label "left/right": Select text for laterality labels, depending on regional practice.

■ L/R:

Display **L** (English: **Left** or German: **Links**) for left label

Display **R** (English: **Right** or German: **Rechts**) for right label

■ G/D:

Display **G** (French: **Gauche**) for left label

Display **D** (French: **Droite**) for right label

■ SIN/DX:

Display **SIN** (Latin: **sinister**) for left label


Display **DX** (Latin: **dexter**) for right label

7.1.7 Deleting images automatically

Image Delete

In the **Image Delete** dialog box you establish the rules for automatic deletion of images in the local database.

Auto delete

Enable auto delete: Activate the automatic image deletion function.

If enabled, you can define the rules for automatic deletion.

- **Hard disc fill level:** Hard disc fill level (depending on the "memory extension" license key) for activation of automatic deletion (60% to 90%)
- **Deletion size:** Amount of hard disk volume to be deleted (10% to 40%)
- **Activation time:** Activation time for automatic deletion (hh:mm:ss)
- **Activation date:** Activation day for automatic deletion (**Daily/Weekly** with selection item for day)

- **Conditions:** Condition(s) for automatic deletion of images (any combination of the following conditions):
 - **Images sent:** Delete all images that have been sent to a DICOM node (with or without confirmation of receipt).
 - **Images filmed:** Delete all printed/filmed images.
 - **Images stored on CD/DVD:** Delete all images stored on CD/DVD.

At least one condition has to be set to activate the automatic delete function.

- **User confirmation:** Obtain user confirmation before starting automatic deletion (show confirmation dialog box).

7.1.8 Password for radiation release

In the **Radiation Release Password** dialog box you change the password required for radiation release.

- 1 Enter the old password.
- 2 Enter a new password (four-digit number).
- 3 Repeat entry of the new password.
- 4 Click OK.

7.2 Examination settings (PEX Editor)

7.2.1 Introduction

The Cios Fusion system provides a series of applications and acquisition program triplets for standard examinations.

The applications are divided into application groups and saved in a database.

The fluoroscopy and acquisition parameters are established in the applications.

All **Fluoro**, **ROAD**, **DR** (Single image), and **SUB** acquisition programs are available as versions with low, medium, and high dose.

Operating mode	Fluoro	ROAD	SUB	DR
Low	x	x	x	x
Medium	x	x	x	x

Operating mode	Fluoro	ROAD	SUB	DR
High	x	x	x	x

An application comprises three acquisition program triplets each (low, medium, high) for each operating mode (fluoro, single image, subtraction, roadmap).

PEX Editor

The **Programs and Exams Editor**, **PEX Editor** for short, is a tool for viewing and editing application groups, applications, and acquisition program triplets. You can also define parameter sets that can be used with the acquisition programs.

The **PEX Editor** may be used by authorized persons only.

PEX database

The PEX Editor administers a database, the **PEX database**, containing the definitions of the application groups, applications, and acquisition program triplets, as well as the parameter sets.

Online/offline

The PEX Editor runs on the Cios Fusion.

Offline editing of the PEX database is also supported. A special version of the PEX Editor runs on a separate Windows PC without the main application.

Data categories

The PEX database supports the storage of the following data types:

■ User data

Your specific applications and application groups.

■ Siemens data

- Standard applications, application groups, parameter sets, and program triplets. These entries are marked with **(d) = default**.
- IQAP applications, application groups, and parameter sets. These entries are marked with **(i) = IQAP**.
- Calibration applications, application groups, and parameter sets. These entries are marked with **(c) = calibration**.
- Service programs for maintenance. These entries are marked with **(m) = maintenance**.

Names

The names of the applications, application groups, parameter sets, and program triplets depend on the language set for the user interface:

- The names of Siemens elements depend on the language set for the user interface. In other words, every element has a different name in each configurable language, but the parameters are the same in each language.
- User specific applications and application groups do not depend on the set language. In other words, there is only one name for all configurable languages.

Storage capacity

The following number of programs can be stored in the database:

- Up to 44 activated application groups, i.e., up to 44 application groups that can be selected on the control panel
- Up to 3000 applications
- Up to 44 applications per application group
- Up to 30 application program triplets per application

Access control

To open the PEX Editor you are asked to enter a password or to select the viewing mode. The customer mode password is site-specific. It can be set up and edited. In service mode, the customer mode password can be reset without knowledge of the current password.

Access modes

There are three modes for accessing the data:

- In **Customer mode**, you have access rights for user data and also for viewing the standard applications and application groups. Access to IQAP, calibration, and service data is not possible.
- In **IQ expert mode**, Siemens image quality experts have access to view and edit all data.
- In **Viewing mode**, you can see all data sets, but you cannot modify them.

Import/Export

You can import a PEX database into the PEX Editor via any local or external drive e.g. local hard drive, a CD/DVD drive, or a USB storage medium. Each import replaces the current database and cannot be undone.

Depending on licensing, several applications/acquisition programs may be inactive after the transfer. Inactive applications/acquisition programs cannot be selected.

Automatic backup

Backups of the PEX database are done automatically. The automatic backup runs in the background after each editing session when the PEX Editor is closed.

Up to ten backups are kept in a folder on the hard drive.

A database stored via automatic backup can be restored with the import function.

7.2.2 Starting the PEX Editor

PREPARATION

1 Change to the **PREPARATION** task card on the control panel if necessary.


2 Press this button.

Manage Tool is displayed on the monitor.

3 Click this icon on the **Settings** subtask card in the **Manage Tool**.

The password dialog box for the **PEX Editor** is displayed.

7 Configuration

The factory-set password for customer mode is: **customer**.

- 4 Select **Customer** mode.
- 5 Enter the correct password.
- 6 Click **Ok** or press **Enter**.


The **PEX Editor** is opened for editing.

– or –

Click **Ok** for **Viewing only**.

The **PEX Editor** is opened for viewing only.

PEX Editor layout

Example

The PEX Editor includes the following task cards:

- **Application Group:** Add or remove applications from application groups.
- **Application:** Establish acquisition program selection and properties within an application.
- **Triplets:** Establish acquisition program triplets for **Fluoro**, **SUB/ROAD**, and **DR** operating modes.
- **Param. Module:** Define and edit parameter sets referenced in the acquisition programs.
- **Configuration:** Define and edit system configuration data, e.g. language configuration, generator type. Functions for backup/restore of the PEX database are also located here.

General work steps in the PEX Editor

- ◆ Check/modify the parameters.
 - (→ Page 224 Editing application groups)
 - (→ Page 225 Editing applications)
 - (→ Page 227 Parameter for fluoroscopy and exposure)

Adopting changed parameters

- ◆ Click **Store**.

Creating a new application group/
application or a new acquisition
program triplet

- 1 Enter a new designation for **Application Group Name** or **Application Name** or **Triplet Name**.
- 2 Click **Store as new**.

7 Configuration


Deleting an application group/
application or an acquisition program
triplet

- 1 Select the element.
- 2 Click **Delete**.
A confirmation dialog box appears.
- 3 Click **Yes** to delete the element.

A triplet used in an application cannot be deleted.
Siemens standard triplets can only be deleted in IQ expert mode.

Rejecting changed parameters

- 1 Click **Close**.

— or —


- Click the Close icon.

A confirmation dialog box appears.

- 2 Click **Ok** to reject changes.

7.2.3 Editing application groups

Example

The **Application Group** task card displays all available applications and application groups. You define application groups and assign applications to the application groups here.

- The applications are listed in the **Application List**.
- The application groups and all assigned applications are listed in the **Application Group List**.

(→ Page 223 *General work steps in the PEX Editor*)

The **Application Group Name** comprises a maximum 15 Unicode characters.

Assigning applications to application groups

- 1 Mark an application in the **Application List**.
- 2 Mark an application group in the **Application Group List**.
- 3 Click **Insert Application**.

– or –

Drag an application from the **Application List** to an application group in the **Application Group List**.

The application is assigned to the application group under the position marked in the **Application Group List**.

The same application can be assigned to multiple application groups.

Removing applications from an application group

- 1 Mark the application in the **Application Group List**.
- 2 Click **Remove Application**.

Establishing the default application group

- ◆ Mark the application in the **Application Group List** and then click **Set as default**.

7.2.4 Editing applications

The screenshot shows the 'Configuration' tab of the 'Application Group' window. The 'Appl. Name' is set to 'Standard (0)'. Below it, there are several rows of settings for different triplet types: Fluoro-Triplet, SUB-Triplet, ROAD-Triplet, DR-Triplet, and 3D-Triplet. Each row has checkboxes for 'Auto Window', 'LH Store', 'Auto-store', 'Storage Rate', 'Max scene length', 'Auto-replay', and 'Number of frames for mask'. The 'Name Module of Contrast Automatic' is set to 'HEG std, DR std, g0.6'. There are also sections for 'Toggle Functions' (EE, MD Threshold, Averaging live, Spatial noise) with 'Low' and 'High' settings, and 'LUT Preset 0' (Center Shift, Width Factor). The 'Flip' section has 'Vert.' and 'Horz.' checkboxes. The 'Magnification' section has a 'Level' slider. At the bottom, there are buttons for 'Store', 'Store as new', 'Delete', 'Previous', and 'Next'.

Example

The **Application** task card shows the assigned acquisition program triplets for the **Fluoro**, **SUB**, **ROAD**, and **DR** operating modes, as well as additional parameters. You can configure the applications and establish the application-specific parameters here.

(→ Page 223 *General work steps in the PEX Editor*)

The **Application Name** comprises a maximum 15 Unicode characters.

Configure triplets for applications

- 1 Select a triplet for each operating mode.
- 2 Activate the **Auto Window** to switch to automatic window value assignment for **SUB/ROAD**.

The brightness and contrast settings of the currently available Peak OP image are used when reviewing the subtraction scene or for subsequent subtraction exposures.
- 3 For automatic saving activate the **Auto-store** option field and enter the **Storage Rate** and **Max scene length**.

The maximum scene length is set in seconds.

The percentage save frequency (save rate) corresponds to the proportion of images to be saved from the live scene; e.g., at a 25% save rate, every 4th image is saved.
- 4 Activate **Autoreplay** if scenes should replay automatically.
- 5 For **SUB/ROAD** establish the **Number of frames for mask**, i.e., the number of images used to calculate the mask.

Defining other parameters

- 1 Select an entry in the **Name Module of Contrast Automatic** list to establish the brightness and contrast.
- 2 Select a multiplication factor for **Toggle function EE**.

The user can toggle the edge enhancement between low, default, and high on the control panel. The default value is set in the triplet. The values for low and high are calculated by a factor.

A factor of 1 means the edge enhancement is medium.
A factor of 0 means the edge enhancement is off.
- 3 Select a multiplication factor for **Toggle function MD**.

The user can toggle the motion detection (k factor) between low, default, and high on the control panel. The default value is set in the triplet. The values for low and high are calculated by a factor.

A factor of 1 means the motion detection is medium.
- 4 Select a **Center Shift** and a **Width Factor** for **LUT Preset 3**.

The values for center shift and width factor can be set for automatic contrast regulation. Correction values for center/width and brightness/contrast of the third LUT can be set.

- 5 Activate **Set Black** if the image area outside the collimated area should be blackened (digital collimation).
- 6 Define the input values for **Rotation** in degrees (0 - 359), **Flip** (vertical and/or horizontal flip) and **Magnification** (Level 0, 1, 2).
- 7 Establish the measurement field with the **Measurement** parameter.

7.2.5 Parameter for fluoroscopy and exposure

Example

The **Triplets** task card enables triplets and parameters to be established for all operating modes on the corresponding cards.

- **Fluoro** card: Parameters for unsubtracted fluoroscopy
- **SUB/ROAD** card: Parameters for subtracted fluoroscopy
- **DR** card: Parameters for unsubtracted single image acquisitions

(→ Page 223 *General work steps in the PEX Editor*)

Parameter

Triplet Name: Name of the fluoroscopy or acquisition program triplet (max. 63 characters).

- **Number of Images:** Number of images from which the LIH (last image hold) is calculated.
- **Number of Images IDL:** Number of pulses used to calculate the individual images in Fluoro and SUB/ROAD modes.
- **Number of Pulses:** Number of pulses to obtain a specific overall dose for a single image acquisition.

Examples: A factor of 1 corresponds to a single acquisition. With a factor of 8, the image is calculated from 8 acquisitions. Possible values: 1, 2, 4, 8, 16, 32

Landmark: Basic setting for anatomical background in subtracted images. Possible values: 0 - 100%

7 Configuration

Parameter modules **Noise Time Filter**, **Noise Spatial Filter**, **Spatial Frequency Response**: Three parameter sets can be selected for each dose level. The parameter sets are defined by experts. (→ Page 228 *Defining parameter sets*)

- **Default:** Establish the default parameter set: **Low** or **Medium**. **High** cannot be selected as the default.

7.2.6 Defining parameter sets

The screenshot shows the 'Param. Module' configuration window with four tabs: 'Spatial Frequency Response Parameter', 'Contrast Automatic Parameter', 'Noise Time Filter Parameter', and 'Noise Spatial Filter Parameter'. Each tab contains various adjustable parameters and buttons for 'Store', 'Store as new', and 'Delete'.

- Spatial Frequency Response Parameter:** Identifier: d009, Name: IQ EE, Harm off. Edge Enhancement: Edge Filter Gain: 2.5, Edge Filter Kernel: 3. Harmonization: Harmonization Gain: 0.0, Harmonization Kernel: 30.0.
- Contrast Automatic Parameter:** Identifier: d010, Name: HEQ std. DR std. g0.6. Peak Specification: MinOP(COF), MaxOP(Jod). NAT: Center: 1,800, Width: 3,500. Histogram: Dynamic Range Limiter: Standard rel, Contrast Limiter: Standard. Auto-gamma: Value: 0.6, Shape: linear. Beta Display: Black, White.
- Noise Time Filter Parameter:** Identifier: d007, Name: IQ k9 MD off. Kmin: 1.0, Kmax: 8.0, Threshold: 300.0. Motion detection: On.
- Noise Spatial Filter Parameter:** Identifier: d004, Name: IQ Noise Filter off. Standard Filter: Level: 1.0. IDEAL Filter: Level: 1.0, Range: 0.0.

Example

Parameter sets can be defined and edited on the **Param. Module** task card. These can then be referenced in application and/or triplet definitions.

Parameter sets may be changed in expert mode only; changed parameter sets cannot be saved in customer mode.

7.2.7 PEX service functions

Configuration parameters

- **Units for:** Unit to be used (cm or inches) for the Collimator Size
- **Password:** (→ Page 229 *Changing the PEX Editor password*)
- **Database Export/Import:** (→ Page 230 *Backing up and restoring the PEX database*)
- **Restore/Import Database Selection:** Database selection for export/import
- **Pex Versions:** Displays the PEX Editor version

Changing the PEX Editor password

- 1 Click the **Configuration** card.
- 2 Click **Change Password**.

The password change dialog is displayed.

- 3 Enter your old password.
- 4 Enter your new password.
- 5 Re-enter your new password.
- 6 Click **Ok**.

Backing up and restoring the PEX database

It is possible to import or export the PEX database via a local or external drive (e.g. local hard drive, CD/DVD drive, USB memory stick).

This process should only be performed by an expert or specialist.

- ◆ Click the **Configuration** card.

Exporting the PEX database


- 1 If necessary, have the required external data medium on hand.

Insert a blank CD/DVD into the drive.

Connect the USB memory to the USB socket.

- 2 Click **Backup / Export**.

A list of all writable devices connected to the PC in use (system or user PC with offline PEX Editor) that can be used for export/backup is displayed.

- 3 Select the destination drive.

- 4 Click **Go** to start export.

The PEX database is copied to the destination drive in a predefined folder (e.g. D:\PEX_backup\). The name of the database copy is the original name extended by creation date and time.

A message window indicates that the backup is complete.

- 5 Click **OK**.

The number of backups on the local drive (C: and D: drives) is limited to 10.

Importing a PEX database By importing you replace the old system database. If an error occurs during transfer, the old database remains valid.

After correct import the old database is permanently deleted.

- 1 If necessary have the required external data medium to which the new database will be stored on hand.

Insert the CD/ DVD into the drive.


Connect the USB memory to the USB socket.

- 2 Select the database you want to replace.

Customer DB, Service DB, or Both for both databases.

- 3 Click **Restore / Import**.

A list of all readable devices connected to the PC in use (system or user PC with offline PEX Editor) that can be used for import/restore is displayed.

- 4 Select the source device.

A list of all restorable databases is shown.

When copying from the local hard drive you have to select one of the two backup folders (manual, automatic backup) to get the list of restorable databases.

- 5 Select the new database or databases.

- 6 Activate **Assure data consistency** if the data are to be checked for errors.

- 7 Click **Go** to start the import.

7 Configuration

8 Maintenance

8.1	Functional and safety checks	234
8.1.1	Daily checks	234
	Prior to the examination	234
	Checking the EMERGENCY STOP function for motor-driven system movements	234
8.1.2	Monthly checks	234
	Functional check of the dose rate control	234
8.1.3	Checks prior to special examinations (e.g. of the open heart and skull)	235
8.1.4	Calibration	235
	Detector calibration	235
	Calibrating the dose measurement chamber	236
8.1.5	Maintenance plan for checking the system	236
8.1.6	Maintenance intervals	237
	General information on maintenance intervals	237
	Batteries	237
8.1.7	Dose and consistency test	238
	Preparation	238
	Dose measurement	238
	Calculation	238
8.2	Service via network connection	239
8.2.1	Remote Service	239
	Software updates via remote access	239
	Setting access rights	239
8.2.2	Application support	241
	Enabling remote access	241
	Rejecting remote access	241
	Entering the session number	241
8.2.3	Transmitting error messages	242
	Transmitting to the factory	242
	Saving to a USB drive	243
8.2.4	Updating virus protection	243
	Behavior during virus attack	243
8.3	Cleaning and Disinfection	244
8.3.1	Cleaning the system parts	244
8.3.2	Cleaning the screen surfaces/TFT displays	245
8.3.3	Disinfection	245
8.4	Protection classes/Protection measures	246
8.5	Active ingredient classes	246

8.1 Functional and safety checks

To ensure that the Cios Fusion is ready for operation and all safety features are functioning properly, you must perform regular functional and safety checks.

8.1.1 Daily checks

Prior to the examination

- 1 Check the power plug. If the power plug is damaged, the Cios Fusion must not be used.
- 2 Check the power cable. If the power cable is damaged, the Cios Fusion must not be used.
- 3 Check the function of the locking brakes of the C-arm system and the monitor trolley as well as the steering of the C-arm system.
- 4 Check the C-arm counterbalance after releasing the brakes.
- 5 Check the loudspeakers for proper functioning. An audible signal must sound while the system starts up.
- 6 Check that all radiation indicators are functioning properly.
- 7 Inspect the flat detector unit housing and the single tank for mechanical damage.
- 8 Check that the scattered radiation grid and the knurled screw are fastened properly.

Checking the EMERGENCY STOP function for motor-driven system movements

- 1 Switch the Cios Fusion on.
- 2 Move the lifting column and press the **EMERGENCY STOP** button at the same time.

The lifting movement is stopped.

A message is displayed on the monitor indicating that EMERGENCY STOP has been actuated.
- 3 Release the **EMERGENCY STOP** button.
- 4 Press one of the buttons for moving the lifting column.

The lifting movement is enabled again.

8.1.2 Monthly checks

Functional check of the dose rate control

You can check the function of the dose rate control without an object using the simple procedure described below. A kV value of ≤ 45 kV must stabilize:

- 1 Open the rectangular and slot collimators to maximum aperture.
- 2 Press the **Tech lock** button and select **Stop**.
- 3 Select 110 kV with the +/- buttons.
- 4 Press the **Tech lock** button again and select **Auto**.

The stop function is canceled; automatic dose rate control is switched on again.

5 Release radiation in fluoroscopy mode.

The tube voltage is reduced to a value ≤ 45 kV.

The monitor image is not overexposed.

8.1.3 Checks prior to special examinations (e.g. of the open heart and skull)

- ◆ Make sure that there is an additional conductive connection between the C-arm system and a point of potential equalization, e.g. the patient table.

8.1.4 Calibration

Detector calibration

The detector calibration shall be done in regular intervals.

CAUTION

Incorrect calibration parameters. Reduced image quality.

Risk of unnecessary radiation exposure due to incorrect calibration parameters!

- ◆ Calibrate the system after 24 months to preserve image quality.

CAUTION

Operator error.

Risk of unnecessary radiation exposure!

- ◆ Follow the calibration procedure instructions.

To ensure optimum image quality, we recommend performing detector calibration at shorter intervals than specified above.

Additionally the imaging system will remind the user to recalibrate the detector in factory-defined intervals.

The detector calibration needs to be done also in case there are new defective pixels visible in the x-ray image. These defective pixels will be displayed e.g. as white dots in the image. In such a case perform the detector calibration as soon as possible.


It will be always possible to work with the system even if the detector calibration warning was confirmed without starting the detector calibration immediately.

- 1 In the **PREPARATION** subtask card, press the **Data system** button on the monitor trolley control panel.

Manage Tool is displayed on the monitor.

8 Maintenance


2 Click the **Service** subtask card on the monitor.

3 Click **Service login**. Alternatively, on the control panel press the **Service login** button on the lower right.

The Service login screen appears.


4 Enter the following password: **admin2** and confirm with **OK**.

The **Tune Up** field appears.

5 Click the **Tune Up** field.

The **Tune Up** screen opens.

Under **User Password** (1) you can enter a separate password or continue with **FD Calibration** (2).

6 Click **FD Calibration** (2).

The FD Calibration screen opens.

7 Follow the menu instructions.

Radiation has to be released during the calibration procedure.

Calibrating the dose measurement chamber

The dose measurement chamber is calibrated on a regular basis as part of the maintenance contract (every 5 years at the latest). If you do not have a maintenance contract, Siemens Customer Service or the manufacturer can calibrate the dose measuring chamber.

8.1.5 Maintenance plan for checking the system

The tests and inspections required by national laws or regulations, such as DHHS regulations or RÖV (constancy tests), are not part of the activities listed in this maintenance plan.

If national laws or regulations specify more frequent checking and/or maintenance, these must be observed.

Maintenance work should be performed by trained technical personnel only. To keep the system in an optimum condition, we recommend that you conclude a maintenance contract. If you have any questions relating to maintenance/ maintenance contract please contact our Siemens Customer Service.

Please observe the relevant information in (→ Page 15 *Maintenance and inspection*).

8.1.6 Maintenance intervals

General information on maintenance intervals

Regular maintenance takes place in 24-month intervals. It contains numerous maintenance steps that have to be performed:

Procedures to be performed	Function	Interval
External visual inspection	Damage, total system	24 months
Mechanical safety	Covers, cable guards, C-arm, C-arm movements, scattered radiation grid/ grid holder, brakes. C-arm system and monitor trolley wheels, lifting column, emergency stop button, warning labels, identification labels, monitors, options (if present, e.g. motorized movements, laser light localizers)	24 months
Electrical safety	Cable and plug, network connection socket, fluoroscopy indicator, acoustic warning signal, radiation release switch, radiation indicator, options (if present, e.g. dose area product and/or air kerma measuring device), total system image quality, checks in accordance with IEC 62353	24 months
Maintenance	System ventilation, cleaning of system	24 months
Operating value inspection	Output and assess Event Log	24 months
Functional inspection	Control functions, monitor indicators of beam path opening/slot collimator positions, options	24 months
Upkeep	Housing	24 months

The stated functions are minimum requirements.

Batteries

Battery replacement intervals (UPS and PC)

Procedures to be performed	Function	Interval
Operating value inspection	UPS lead gel batteries PC bios battery	Every 48 months

The stated functions are minimum requirements

8.1.7 Dose and consistency test

Carry out dose area product testing after performing the following procedure:


Preparation

- 1 Open the collimator completely (full format).
- 2 Select the "Fluoroscopy" operating mode.
- 3 Center the lead ruler on the flat detector.
- 4 Release radiation, setting the rectangular collimator to approx. 15 cm edge length.
- 5 Keep the field size (shown on the lead ruler) constant.

The rectangular collimator may no longer be adjusted.

- 6 Remove the lead ruler.
- 7 Center the dose measurement device on the flat detector (see (1) in graphic).

(1) Dose measurement device

Dose measurement

- 1 Release radiation for a short period and then actuate the **Teck lock** button and set 70 kV.
- 2 Record the dose area product - value "A" (μGycm^2) - shown on the monitor.
- 3 Set the measured dose area product to "0" on the measuring device.
- 4 Release radiation for approx. 10 s and record the measured dose (cGy).
- 5 Record the dose area product - value "B" - shown on the monitor.
- 6 Subtract value "A" from value "B".

Calculation

Measured dose area product = measured dose (cGy) x 225 cm^2

- ◆ Compare the measured value with the dose area product displayed on the monitor.

With a correct adjustment, the displayed and calculated values agree within a maximum deviation of +/- 25%. If the deviation exceeds this percentage, please inform service.

8.2 Service via network connection

The Cios Fusion must be connected to a network (DSL $\geq$ 1 Mbit/s) to enable remote access service measures.

8.2.1 Remote Service

If the rights are granted accordingly, Siemens Customer Service can access the Cios Fusion for maintenance and installation activities via an active network connection.

Software updates via remote access

CAUTION

Software installation interferes with system functions.

Interruption of patient examination and loss of patient data.

- ◆ Do not start the installation while working with patients.

CAUTION

Failure installing the update package.

A failed software installation means that the system is in an undefined state.

- ◆ Stop using the system and notify the UPTIME Service Center.

CAUTION

No restart after service session.

Changes to the service environment may not go into effect without a system restart.

- ◆ Restart the system after a service session.

Setting access rights

In the **Remote Service** window you grant service personnel access rights so that they can perform maintenance in a remote session.

The patient must be closed; otherwise, there could be problems accessing the control panels.

- 1 In the **PREPARATION** subtask card, press the **Data system** button on the monitor trolley control panel.


Manage Tool is displayed on the monitor.

- 2 Click the **Settings** subtask card on the monitor.

- 3 Double-click the Telephone icon.

The **Remote Service** window appears.

- 4 Click one of the following options to define the type of access:

Full access: Service receives full access rights to your system.

Limited access: Service receives limited access, that is, they can see error messages and reset passwords.

No access (default setting): Service does not receive access to your system.

- 5 Click **OK** to apply all settings and exit the window.


- 6 Press this button again to exit the **Manage Tool**.

8.2.2 Application support

If you require online support for your system, you can grant access to the Cios Fusion as follows:

- 1 Right-click the **i icon** (1) on the monitor.

The **Headquarter Support Applications** screen opens:

- 2 Click **Remote assistance**.

The **Remote Desktop Assistance** window opens.

Enabling remote access

- ◆ To enable remote access to the system click **I agree**.

The dialog window for entering the session number opens.

Rejecting remote access


- ◆ To reject remote access to the system click **I do not agree**.

The request is canceled.


Entering the session number

- 1 Enter the session number provided to you by Service.

8 Maintenance

2 Click **Login**.

- 3 Click **OK** to accept the request to take over your system remotely.
The requesting technician receives the right to control the user interface.
- 4 Click **Viewing Only**, to deactivate remote system control if it is currently active.
System control by the moderator or other participant is cancelled. System control is returned to you.
- 5 Click **Disconnect** to disconnect the established connection.
The remote support connection icon disappears.
- 6 Click **Cancel** to close the window.


8.2.3 Transmitting error messages

You can transmit error messages directly to the factory or save them to a USB drive.

Transmitting to the factory

1 Click the **i** icon.

The **Headquarter Support Applications** screen opens.

2 Under **Send to factory** mark the appropriate check box.

3 You can also enter a comment by marking the appropriate check box and entering text.

4 Confirm with **OK**.

The message and/or log files are sent.


If the system is not currently connected to the network, the files are placed in a transfer folder on the system. The files are sent the next time a network connection is established.

Saving to a USB drive

1 Insert the USB drive into a free USB slot on the monitor trolley.

2 Click the corresponding box and confirm with **OK**.

The log files will be saved to the USB drive.

8.2.4 Updating virus protection

A virus scanner is part of the software installed on your system.

CAUTION

New viruses not detected.

System malfunction as the result of a virus attack.


- ◆ Acknowledge the virus pattern file and scan engine updates during booting. Only available with remote handling of updates.
- ◆ Otherwise, contact Siemens Customer Service.

The Cios Fusion should be operated only in a virus-protected hospital network.

Behavior during virus attack

A message window is displayed if a virus is detected.

In this case, you must contact Siemens Customer Service.
Do not try to correct the problem yourself.

8.3 Cleaning and Disinfection

Before cleaning or disinfecting the Cios Fusion, the system must be disconnected from the power supply and switched off.

CAUTION

Inadvertent actuation of the foot pedal.

Risk of radiation!

- ◆ Switch off the system before cleaning.

Before each examination, clean all parts that could come into contact with the patient or bodily fluids to prevent contamination of the Cios Fusion. Also make sure that any dust deposits on the housing of the image intensifier/flat detector, the single tank, and the C-arm are removed prior to every treatment. Please follow the cleaning instructions as described below.

CAUTION

Improper cleaning.

Risk of infection!

- ◆ After dealing with patients with highly infectious diseases, clean the unit wheels, in addition to regular cleaning.

8.3.1 Cleaning the system parts

- 1 Clean the system parts with a damp cloth.

- 2 For moistening, use water or a lukewarm, diluted aqueous solution consisting of water and a household cleaning agent.

Cleaning agents and disinfectants containing chlorine derivatives can discolor cables, but will not negatively impact system functions.

Never immerse system parts (excepting the footswitch) in liquid or autoclave system parts!
Observe the respective protection and active ingredient classes in this regard.

8.3.2 Cleaning the screen surfaces/TFT displays

The screen surfaces can be cleaned at any time.

- 1 Clean the monitor screen with a cotton cloth dampened with water.
- 2 Immediately dry off the monitor screen with a soft cotton cloth.
- 3 Wipe off contrast agent spots as soon as possible.

8.3.3 Disinfection

For the disinfection of surfaces we recommend liquid solutions of common surface disinfectants based on aldehyde and/or amphoteric surfactants, e.g. Tensodur 103, Kohrsolin, Cidex.

Certain substituted phenol-based or chlorine-splitting disinfectants can corrode materials and are therefore not recommended.

CAUTION

Use of unsuitable cleaning agents.

Risk of contamination.

- ◆ Only use the recommended cleaning agents and disinfectants.

CAUTION

Use of abrasive cleaning agents, liquids, or sprays.

Risk of electrical shock or damage to the system.

- ◆ Only use the recommended cleaning agents and disinfectants.
- ◆ Do not allow any cleaning fluids to enter the openings of the system (e.g. air openings, gap between covers).
- ◆ Follow the cleaning and disinfection instructions.

- ◆ Please also observe the instructions for use of the disinfectant.

Some substances contained in disinfectants are known to be hazardous to health. The concentration of such substances in the air must not exceed the legally defined limit. We recommend that you follow the manufacturers' usage instructions for these products (active ingredient classes).

8.4 Protection classes/Protection measures

System components	Protection class	Protection measures
X-ray tube assembly (Single-tank)	IPX0 - No protection	protect with sterile cover
Flat detector	IPX0 - No protection	protect with sterile cover
C-arm	IPX0 - No protection	protect with sterile cover
Remote control unit (optional)	IP65 dust tight, protected against water jets	protect with sterile cover, separate from chassis before cleaning
Standard footswitch	IPX8 protected against immersion beyond 1 m	protect with sterile cover, separate from chassis before cleaning
Multifunctional footswitch	IPX8 protected against immersion beyond 1 m	protect with cover, separate from chassis before cleaning
Hand switch	IP54 dust protected, protected against splashing of water	separate from chassis before cleaning

8.5 Active ingredient classes

The following active ingredient classes may be used to clean and disinfect the Cios Fusion.

Active ingredient class	Aldehydes, alkylamines, quaternary compounds, guanidine derivatives, peroxide compounds, organic acids, chlorine derivatives, alcohol, benzene, common household dish-washing liquids, detergents
-------------------------	---

9 Technical Specifications


9.1	Curves and diagrams	249
9.1.1	Heating and cooling curves of the single tank	249
9.1.2	X-ray tube emission curve	249
9.1.3	Dose rate at the flat detector input	249
	Dose rate value deviation	250
	Setting the dose rate value	250
9.1.4	Fluoro curves	250
	Systems with 2.3 kW generator	250
9.1.5	Dosimetric information	251
	Air kerma strength	251
	Dosimetric information in accordance to IEC 60601-2-43:2010	259
9.2	Notes on electromagnetic compatibility (EMC)	261
9.2.1	Guidelines and manufacturer's declaration – Electromagnetic emissions	261
9.2.2	Guidelines and manufacturer's declaration – Electromagnetic interference immunity	262
9.2.3	Guidelines and manufacturer's declaration – Electromagnetic interference immunity	263
9.2.4	Recommended separation distances between portable and mobile RF telecommunications equipment and the system	264
9.3	Original equipment manufacturer	265
9.3.1	Introduction	265
9.3.2	For Open Source Software (OSS) only	265
9.4	System data	266
9.4.1	Entire system	266
	General data	266
	Weight	267
	Classification	267
	Current/voltage values	267
9.4.2	Generation of radiation	267
	X-ray generator	267
	X-ray tube assembly (Single-tank)	268
9.4.3	Image receptor	269
	Flat detector	269
	Scattered radiation grid on flat detector (detachable)	271
9.4.4	Unit components	271
	C-arm	271
	Collimator system	272
	Imaging chain	272
	Monitors	272
9.4.5	Equipment options	273
	Dose measuring chamber (DAP meter)	273
	Automatic dose rate control (ADR)	274
	Detector laser light localizer (optional)	274
	Single-tank laser light localizer (optional)	274
	HD video manager	274
	DVI video splitter	274
	SCALANCE W 721 WLAN access point	275
9.5	List of Standards	275

9 Technical Specifications

9.6	Labels	276
9.6.1	C-arm system	277
	1. Drawer on C-arm travel frame	277
	2. Base system (chassis)	280
	3. C-arm system control unit and remote control unit	281
	4. Flat detector, grid, frame, laser light localizer	281
	5. Single tank	282
	6. Single tank laser and dose measuring chamber (optional)	283
9.6.2	Monitor trolley	284


9.1 Curves and diagrams

9.1.1 Heating and cooling curves of the single tank

- (1) blue curve: Cooling curve
- (2) orange curve: Heating curve (80 W)
- (3) green curve: Heating curve (330 W)
- (4) brown curve: Heating curve (600 W)

9.1.2 X-ray tube emission curve

9.1.3 Dose rate at the flat detector input

The dose rate is set at the factory; it is measured behind the scattered radiation grid at the flat detector input. Refer to the acceptance protocol §16 RÖV for the system-specific dose values (Federal Republic of Germany only).

The dose rate is set in the kV range between 70 and 80 kV using a technical phantom.

Dose rate value deviation

Depending on the examined object, different fluoroscopy data (kV, mA) arise. Because flat detector sensitivity depends on beam quality (kV), the same gray values in the image may correspond to different dose rate values at the detector input.

When examining a patient in fluoroscopy mode, additional scattered radiation values are produced in comparison to the phantom values, affecting the dose rate at the flat detector input.

Setting the dose rate value

If desired, the preferred position for the dose rate can be reprogrammed.

9.1.4 Fluoro curves

The following fluoro curve properties are available in the system for clinical use.

The following base values apply for all systems and fluoro curves:

- Minimum kV = 40
- Maximum kV = 110

The currents (mA) shown in the following tables indicate effective values.

Systems with 2.3 kW generator

Characteristic curves for fluoroscopy

The maximum average power for normal fluoroscopy characteristic curves is 600 W. Characteristic curves with the HL add-on have a maximum average power of 1000 W (limited to 30 s radiation).

Characteristic curve	Min mA (pulse)	Max mA (pulse)	Example for the same object at 15 fps ¹⁾		
			kV	Ø mA (calculated)	Ø power (W)
FL-CARD	3.0	23.7	80	4.0	320
FL-CARD HL	3.0	23.7	80	4.0	320
FL-HC	3.0	23.7	69	7.9	550
FL-HC HL	3.0	23.7	68	8.4	574
FL-IOD	3.0	23.7	67	8.9	599
FL-IOD HL	3.0	23.7	66	9.5	625
FL-LD	3.0	23.7	92	2.2	205
FL-S	3.0	23.7	85	3.0	254
FL-S HL	3.0	23.7	82	3.6	292

1) with the same image receptor dose

Characteristic curves for single image

Characteristic curve	Min mA (pulse)	Max mA (pulse)	Example for the same object ²⁾	
			kV	mAs
DR-CARD	3.0	23.7	80	0.266
DR-HC	3.0	23.7	68	0.562
DR-IOD	3.0	23.7	66	0.630
DR-LD	3.0	23.7	86	0.188
DR-S	3.0	23.7	66	0.631

9.1.5 Dosimetric information

Air kerma strength

The skin dose values were measured at a distance of 30 cm (patient entry reference point) from the image receptor input with a 20 cm PMMA phantom (equivalent to a typical patient) in accordance with IEC 60601-2-43:2010, 203.5.2.4.5.101.

Radioscopy, Detector size 20x20

Program ³⁾ with scattered radiation grid	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate display [mGy/min]
Intervention, Card Interv. Zoom 0, 20 x 20	low	80	9.4	7.1	15	3.4
	medium	80	2.9	x	cont.	10.7
	high	80	4.5	x	cont.	16.1
Intervention, Card Interv. Zoom 1, 15 x 15	low	80	15.0	7.1	15	5.5
	medium	80	4.0	x	cont.	14.5
	high	80	5.8	x	cont.	20.7
Intervention, Card Interv. Zoom 2, 10 x 10	low	80	20.0	7.1	15	7.5
	medium	80	6.2	x	cont.	22.6
	high	80	9.4	x	cont.	34.2
Intervention, Card Gastro Zoom 0, 20 x 20	low	91	8.8	7.1	15	3.4
	medium	91	1.6	x	cont.	7.6
	high	94	2.1	x	cont.	10.9

2) with the same image receptor dose

3) Only valid as delivered from factory. Technical changes reserved.

9 Technical Specifications

Program ³⁾ with scattered radiation grid	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate display [mGy/min]
Intervention, Card Vascular	low	63	24.0	12.1	7.5	4.3
	medium	66	24.0	13.5	15.0	10.8
	high	66	24.0	19.4	15.0	16.0
Cardio Card Standard Zoom 0, 20 x 20, cont.	low	80	2.0	x	cont.	7.2
	medium	80	2.9	x	cont.	10.7
	high	80	4.5	x	cont.	16.7
Gastro Standard Zoom 0, 20 x 20, cont.	low	88	1.3	x	cont.	5.7
	medium	91	1.6	x	cont.	7.7
	high	94	2.1	x	cont.	10.9
Vascular, Vasc Standard Zoom 0, 20 x 20	low	64	24.0	10.4	15	7.7
	medium	66	24.0	13.3	15	10.6
	high	67	24.0	17.3	15	14.8
Urology, Uro Standard Zoom 0, 20 x 20	low	64	24.0	10.4	15	7.7
	medium	66	24.0	13.3	15	10.6
	high	67	24.0	17.3	15	14.8
Ortho/Trauma Ortho Standard Zoom 0, 20 x 20	low	64	24.0	10.4	15	7.7
	medium	66	24.0	13.3	15	10.6
	high	67	24.0	17.3	15	14.8
General Standard Zoom 0, 20 x 20	low	64	24.0	10.4	15	7.7
	medium	66	24.0	13.3	15	10.6
	high	67	24.0	17.3	15	14.8
General Standard HL Zoom 0, 20 x 20	low	80	9.4	7.1	15	3.4
	medium	80	18.0	7.1	15	6.5
	high	80	24.0	7.1	15	17.8
General, Standard HL, high dose, cont.	high	110	8.9	x	cont.	61.8

3) Only valid as delivered from factory. Technical changes reserved.

Program ³⁾ with scattered radiation grid	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate display [mGy/min]
General Pediatrics	low	69	6.3	7.1	15	1.5
w/o scattered radiation grid	medium	73	6.8	7.1	15	1.9
	high	78	7.4	7.1	15	2.5

Radioscopy, Detector size 30x30

Program ⁴⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate display [mGy/min]
Intervention, Card Interv. Zoom 0, 30 x 30	low	77	7.2	7.1	15	2.0
	medium	80	1.9	x	cont.	5.8
	high	80	2.8	x	cont.	8.8
Intervention, Card Zoom 1, 20 x 20	low	80	9.4	7.1	15	3.1
	medium	80	3.0	x	cont.	9.7
	high	80	4.5	x	cont.	15.8
Intervention, Card Zoom 2, 15 x 15	low	80	13.0	7.1	15	4.1
	medium	80	4.0	x	cont.	13.2
	high	80	6.0	x	cont.	19.4
Intervention, Card Gastro Zoom 0, 30 x 30	low	76	7.10	7.1	15	1.9
	medium	88	1.2	x	cont.	4.5
	high	90	1.5	x	cont.	6.0
Intervention, Card Vascular Zoom 0, 30 x 30	low	60	24.0	7.6	7.5	2.0
	medium	63	24.0	9.2	15.0	5.7
	high	64	24.0	13.0	15.0	8.3
Cardio Card Standard Zoom 0, 30 x 30, cont.	low	80	1.3	x	cont.	3.9
	medium	80	1.9	x	cont.	5.9
	high	80	2.8	x	cont.	8.9

3) Only valid as delivered from factory. Technical changes reserved.

4) Only valid as delivered from factory. Technical changes reserved.

9 Technical Specifications

Program ⁴⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate display [mGy/min]
Gastro Standard Zoom 0, 30 x 30, cont.	low	85	0.9	x	cont.	3.2
	medium	88	1.2	x	cont.	4.5
	high	90	1.5	x	cont.	6.0
Vascular, Vasc Standard Zoom 0, 30 x 30	low	62	24.0	7.1	15	4.0
	medium	63	24.0	9.3	15	5.8
	high	65	24.0	12.1	15	8.1
Urology, Uro Standard Zoom 0, 30 x 30	low	62	24.0	7.1	15	4.0
	medium	63	24.0	9.3	15	5.8
	high	65	24.0	12.1	15	8.1
Ortho/Trauma Ortho Standard	low	62	24.0	7.1	15	4.0
	medium	63	24.0	9.2	15	5.8
	high	65	24.0	12.1	15	8.1
General Standard Zoom 0, 30 x 30	low	62	24.0	7.1	15	4.0
	medium	63	24.0	9.3	15	5.8
	high	65	24.0	12.1	15	8.1
General Standard HL Zoom 0, 30 x 30	low	80	7.5	7.1	15	2.3
	medium	80	14.0	7.1	15	4.5
	high	80	24.0	11.2	15	12.4
General, Standard HL, high dose, cont.	high	110	8.9	x	cont.	54.5
General Pediatrics w/o scattered radiation grid	low	64	5.8	7.1	15	0.9
	medium	68	6.3	7.1	15	1.2
	high	72	6.7	7.1	15	1.6

Other pulse rates and zoom factors Dose values for other pulse rates and zoom factors may be roughly estimated with the aid of the following scaling factors:

4) Only valid as delivered from factory. Technical changes reserved.

Pulse rate [1/s]	Scaling Factors for Air Kerma/Air Kerma Rate ⁵⁾ related to 15 p/s	
0.5	0.03	
1	0.07	
2	0.13	
3	0.2	
5	0.33	
7.5	0.50	
10	0.67	
15	1.00	
30	2.00	

Zoom Steps [Magnify]	Scaling Factors for Air Kerma/Air Kerma Rate ⁶⁾	
	20x20 Detector	30x30 Detector
0	1	1
1	1.33	1.5
2	2	2

Radiography, Detector size 20x20

Program ⁷⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Charge [mAs]	Dose [mGy]
Intervention, Card Zoom 0, 20 x 20	low	80	1.50	0.10
	medium	99	1.78	0.20
	high	110	1.59	0.20
Intervention, Card Zoom 1, 15x15	low	80	1.88	0.10
	medium	110	1.59	0.20
	high	110	1.59	0.20

5) Only valid as delivered from factory. Technical changes reserved.

6) Only valid as delivered from factory. Technical changes reserved.

7) Only valid as delivered from factory. Technical changes reserved.

9 Technical Specifications

Program ⁷⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Charge [mAs]	Dose [mGy]
Intervention, Card Zoom 2, 10x10	low	80	1.0	< 0.10
	medium	110	1.59	0.20
	high	110	1.59	0.20
Intervention, Card Gastro Zoom 0, 20x20	low	93	0.75	< 0.10
	medium	99	1.78	0.20
	high	110	1.59	0.20
Intervention, Card Vascular Zoom 0, 20x20	low	76	2.41	0.10
	medium	108	1.63	0.20
	high	110	1.59	0.20
Card, Card Standard Zoom 0, 20x20, cont.	low	80	1.50	0.10
	medium	99	1.78	0.20
	high	110	1.59	0.20
Gastro Standard	low	93	0.75	< 0.10
	medium	99	1.78	0.20
	high	110	1.59	0.20
Vasc. Vascular Standard Zoom 0, 20x20	low	74	2.47	0.10
	medium	105	1.68	0.20
	high	110	1.59	0.20
Urology, Uro Standard Zoom 0, 20x20	low	74	2.47	0.10
	medium	105	1.68	0.20
	high	110	1.59	0.20
Ortho/Trauma, Ortho Standard Zoom 0, 20x20	low	74	2.47	0.10
	medium	105	1.68	0.20
	high	110	1.59	0.20
General Standard Zoom 0, 20x20	low	74	2.47	0.10
	medium	105	1.68	0.20
	high	110	1.59	0.20

7) Only valid as delivered from factory. Technical changes reserved.

Program ⁷⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Charge [mAs]	Dose [mGy]
General Standard HL Zoom 0, 20x20	low	80	1.50	< 0.10
	medium	99	1.78	0.20
	high	110	1.59	0.20
General, Standard HL, high dose, cont.	high dose	110	1.59	0.20
General Pediatrics w/o scattered radiation grid	low	88	0.31	< 0.10
	medium	93	0.71	< 0.10
	high	96	1.24	< 0.10

Radiography, Detector size 30x30

Program ⁸⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Charge [mAs]	Dose [mGy]
Intervention, Card Zoom 0, 30 x 30	low	80	0.89	< 0.10
	medium	86	2.06	0.20
	high	108	1.63	0.20
Intervention, Card Zoom 1, 20x20	low	80	1.00	< 0.10
	medium	99	1.78	0.20
	high	110	1.59	0.20
Intervention, Card Zoom 2, 15x15	low	80	0.89	< 0.10
	medium	110	1.59	0.20
	high	110	1.59	0.20
Intervention, Card Gastro Zoom 0, 30x30	low	90	0.47	< 0.10
	medium	96	1.19	0.10
	high	105	1.68	0.20
Intervention, Card Vascu- lar Zoom 0, 30x30	low	72	1.56	< 0.10
	medium	87	2.08	0.20
	high	110	1.59	0.20

7) Only valid as delivered from factory. Technical changes reserved.

8) Only valid as delivered from factory. Technical changes reserved.

9 Technical Specifications

Program ⁸⁾ with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Charge [mAs]	Dose [mGy]
Card, Card Standard Zoom 0, 30x30, cont.	low	80	0.89	< 0.10
	medium	86	2.06	0.20
	high	108	1.63	0.20
Gastro Standard Zoom 0, 30x30, cont.	low	90	0.47	< 0.10
	medium	96	1.19	0.10
	high	105	1.68	0.20
Vasc. Vascular Standard Zoom 0, 30x30	low	72	1.56	< 0.10
	medium	87	2.08	0.20
	high	110	1.59	0.20
Urology, Uro Standard Zoom 0, 30x30	low	72	1.56	< 0.10
	medium	87	2.08	0.20
	high	110	1.59	0.20
Ortho/Trauma, Ortho Standard Zoom 0, 30x30	low	72	1.56	< 0.10
	medium	87	2.08	0.20
	high	110	1.59	0.20
General Standard Zoom 0, 30x30	low	72	1.56	< 0.10
	medium	87	2.08	0.20
	high	110	1.59	0.20
General Standard HL Zoom 0, 30x30	low	80	0.95	< 0.10
	medium	99	1.93	0.20
	high	110	1.59	0.20
General, Standard HL, high dose, cont.	high dose	110	1.59	0.20
General Pediatrics w/o scattered radiation grid	low	86	0.20	< 0.10
	medium	90	0.49	< 0.10
	high	94	0.86	0.10

8) Only valid as delivered from factory. Technical changes reserved.

Dosimetric information in accordance to IEC 60601-2-43:2010

- Measuring conditions**
- Specifications for interventional systems acc. to EN/IEC 60601-2-43
 - Reference air kerma strength and ranges of air kerma strength according to EN/IEC 60601-2-43

Isokerma curves


All following curves were recorded as follows:

beam projection horizontal and vertical (tube unit above)

- SID of 102.5 cm
- 25 cm x 25 cm x 25 cm PMMA Phantom centered/aligned central beam axis
- Entrance Surface of the Phantom is at the PATIENT ENTRANCE REFERENCE POINT
- Radiation field size 10 cm x 10 cm at entrance of PMMA Phantom
- Isokerma curves 100 cm and 150 cm above the floor
- Scattered radiation of isokerma curves in $\mu\text{Gy}/(\text{Gy cm}^2)$ at 110 kV / 5.3 mA continuous according to IEC 60601-2-43:2010


The values in $\mu\text{Gy}/(\mu\text{Gy m}^2)$ according to IEC 60601-2-43:2000 can be calculated by a factor of 0.01.

Cios Fusion, 100 cm above floor,
horizontal, $\mu\text{Gy}/(\text{Gy cm}^2)$


9 Technical Specifications


Cios Fusion, 100 cm above floor,
vertical, $\mu\text{Gy}/(\text{Gy cm}^2)$

Cios Fusion, 150 cm above floor,
horizontal, $\mu\text{Gy}/(\text{Gy cm}^2)$

Cios Fusion, 150 cm above floor,
vertical, $\mu\text{Gy}/(\text{Gy cm}^2)$

9.2 Notes on electromagnetic compatibility (EMC)

Medical electrical equipment is subject to special precautions regarding EMC. These systems must be installed and put into service according to the EMC information provided in the accompanying documents.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Permanently connected equipment or system cabling that cannot be removed by the user is not listed. This cabling is part of the system and was considered in all EMC measurements. Without this cabling the equipment or system would not function.

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emission or decreased immunity of the equipment or system.

9.2.1 Guidelines and manufacturer's declaration – Electromagnetic emissions

The system is intended for use in the electromagnetic environment as specified below. The customer or the user of the system should ensure that it is operated in such an environment.

9 Technical Specifications

Emissions test	Compliance	Electromagnetic environment – Guidelines
RF emissions according to CISPR 11	Group 1	The system uses RF energy exclusively for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions according to CISPR 11	Class A	The system is not suitable for use in domestic establishments and establishments directly connected to a public power supply network that supplies buildings used for domestic purposes.
Harmonic emissions according to IEC 61000-3-2	Not applicable	The system is a professionally used unit with a total rated power of more than 1 kW. There are no limit values for this rated power.
Emission of voltage fluctuations/flicker according to IEC 61000-3-3	Complies	

The use of this device or system can cause RF interference or interfere with the operation of nearby equipment. In such cases employ suitable remedies such as shielding or rearranging the equipment.

The device or system should not be operated next to or on/under other devices. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

9.2.2 Guidelines and manufacturer's declaration – Electromagnetic interference immunity

The system is intended for use in the electromagnetic environment as specified below. The customer or the user of the system should ensure that it is operated in such an environment.

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment – Guidelines
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic materials, the relative humidity must be at least 30%.
Electrical fast transient/bursts according to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	The line power quality should correspond to that of a typical commercial or hospital environment.

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment – Guidelines
Surges according to IEC 61000-4-5	± 1 kV differential-mode voltage ± 2 kV common-mode voltage	± 1 kV differential-mode voltage ± 2 kV common-mode voltage	The line power quality should correspond to that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the power supply voltage according to IEC 61000-4-11	< 5% U_T * (> 95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	< 5% U_T * (> 95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	The line power quality should correspond to that of a typical commercial or hospital environment. If the user of the system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptible power supply.
Magnetic field at supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a commercial or hospital environment.

* Remark: U_T is the AC line power voltage before application of the test level.

9.2.3 Guidelines and manufacturer's declaration – Electromagnetic interference immunity

The system is intended for use in the electromagnetic environment as specified below. The customer or the user of the system should ensure that it is operated in such an environment.

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment – Guidelines
			Portable and mobile RF telecommunications equipment should be used no closer to any part of the system (including cables) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF disturbances according to IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	10 Vrms	$d = 1.2 \sqrt{P}$

9 Technical Specifications

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment – Guidelines
Radiated RF disturbances according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz
			$d = 2.3\sqrt{P}$ for 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by electromagnetic site surveya, should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 
Remark 1: At 80 MHz and 800 MHz, the safe distance for the higher frequency range applies.			
Remark 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To be able to assess the electromagnetic environment with regard to the stationary transmitters, an electromagnetic site study should be considered. If the measured field strength in the location in which the device or system is used exceeds the applicable RF compliance level above, the device or system should be observed to verify normal operation. Should unusual performance features be observed, additional measures (such as change in orientation or change of site of the system) may be necessary.			
b Above the frequency range of 150 kHz to 80 MHz the field strength should be less than 3 V/m.			

9.2.4 Recommended separation distances between portable and mobile RF telecommunications equipment and the system

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the operator of the system can contribute to avoiding electromagnetic interference. For this purpose he should comply with the minimum distance stated below (depending on the maximum output nominal power of the communication devices) between portable and mobile RF telecommunication devices (transmitters) and the system.

Rated maximum output power of transmitter [W]	Safe distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed in the table above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Remark 1: At 80 MHz and 800 MHz, the safe distance for the higher frequency range applies.			
Remark 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.			

9.3 Original equipment manufacturer

Original Equipment Manufacturer (OEM), Open Source Software (OSS) from third-party providers in Siemens Healthcare products.

9.3.1 Introduction

The product can contain OEM and OS software licensed and developed by third-party providers (depending on the system configuration and options).

These OEM and/or OS files are protected by copyright. Your right to use OEM and/or OS software beyond the mere execution of the Siemens program is governed by the relevant terms of the OEM and/or OS software license.

Complying with these licensing conditions authorizes you to use the OEM and/or OS software as provided in the respective license. In case of conflicts between the licensing conditions of Siemens and the OEM and/or OS software license conditions, the OEM and/or OS software conditions govern the OEM and/or OS portions of the software.

The OEM and/or OS software licenses (license conditions) are available on the CDs/DVDs (open source/legal concept) delivered together with the system.

9.3.2 For Open Source Software (OSS) only

Open source software is licensed free of charge (meaning that exercising the licensing rights is free of charge, while fees may be charged to recover the costs borne by SIEMENS). If programs contained in this product are licensed under GNU General Public License (GPL), GNU Lesser General Public License (LGPL), or another open source license that requires source code, the corresponding code is delivered together with the device on CD/DVD.

Warranty regarding the further use of software from third-party providers:

SIEMENS makes no warranty regarding the OEM and/or OS software programs in this device, when such programs are used in a manner other than the program execution intended by SIEMENS. The licenses mentioned previously define the warranty, if any, from the authors or licensors of the OEM and/or OS software. SIEMENS specifically disclaims any warranties for defects caused by altering any OEM and/or OS software program or the product's configuration. If the OEM and/or OS software violates the intellectual property of a third party, you have no warranty claim against SIEMENS.

Technical support is available only for software that has not been modified.

9.4 System data

9.4.1 Entire system

General data

Power requirements	100 V, 110 V, 120 V, 127 V, 200 V, 220 V, 230 V, 240 V $\pm 10\%$; 50/60 Hz (± 1 Hz) An on-site 16 A/characteristic C overcurrent release (trip breaker) is recommended by DIN VDE 0100-710 for supply network voltages of 230 V/240 V. Country-specific requirements apply for the fuse values to be used for other supply network voltages.
Unit fuse protection (internal)	100 V to 127 V 20 A slow-blow fuse 200 V to 240 V 15 A slow-blow fuse
Nominal rating	20 A to 127 V~, 15 A from 200 V~ corresponds to the nominal value of the slow-blow fuse in the power input of the product
Internal line impedance	Ri max. 0.3 Ohm for 100 - 127 V~ Ri max. 0.8 Ohm for 200 - 240 V~
Power consumption	Maximum 2.6 kW Standby 475 W (for 230 V)
Continuous	14.5 A (100 V) / 7 A (230 V)
Short term (max. 30 s 1000 W Fluoro)	25 A (100 V) / 12 A (230 V)
Short term (max. ~ 50 ms DR single shot)	29 A (100 V) / 14 A (230 V)
Ambient conditions for operation	Temperature range: +15 °C to +35 °C Rel. humidity: 15% to 75%, non-condensing Barometric pressure: 700 hPa to 1060 hPa

Ambient conditions for transport and storage	Transport temperature range: -25 °C to +45 °C Storage temperature range: -25 °C - +37 °C Rel. humidity: 10% to 95%, non-condensing Barometric pressure: 700 hPa to 1060 hPa
--	---

For ambient conditions that fall within the specified values, no waiting time is required for radiation operation.

For ambient conditions that fall outside the specified values, a waiting time of one to twelve hours must be factored in for radiation operation, depending on the location of the system.

Weight

Monitor trolley (including 2 monitors, UPS)	230 kg
C-arm system without accessories	272 kg

Classification

Protection against electric shock	Class 1, no applied part according to IEC 60601-1
Protection against ingress of fluids	IPX0 (not protected) according to IEC 60529
Operating mode	Pulsed and continuous radiology

Current/voltage values

Voltage (V AC)	Long-term current consumption 50 Hz/60 Hz (A)	Short-term current consumption 50 Hz/60 Hz (A)
100	18	29
110	17	27
120	16	25
127	16	24
200	10	16
220	9	14
230	9	14
240	9	14

9.4.2 Generation of radiation

X-ray generator

The values in the table provide the technical thresholds of X-ray generation. Whether these threshold values can be reached during clinical operation depends on the set configuration values.

9 Technical Specifications

High-frequency generator	2.3 kW (non-IEC)
Rated power (IEC 60613)	1.4 kW (100 kV, 14 mA, 100 ms)
Single image pulse power (max)	1.9 kW (80 kV, 23.7 mA, 33.5 ms)
Average power in series operation (max)	1000 W (71 kV, 14.1 mA, 33.3 ms, 30 P/s)
Inverter control frequency	15 kHz to 30 kHz
kV range	40 kV to 110 kV
Single image	0.2 mA to 25 mA
	0.021 mAs to 4 mAs
	7 ms to 160 ms
Continuous fluoroscopy	0.2 mA to 15 mA
Fluoroscopy	3 mA to 25 mA
	0.5 P/s to 15 P/s (30 P/s for continuous fluoroscopy)
	7 ms to 40 ms pulse width
Single tank heat storage capacity	888,000 J/1,200,000 HU

Data according to IEC 60601-2-54:

Radiography	
Rated tube voltage	110 kV
Highest X-ray tube current	25 mA
Highest power (nominal power)	100 kV, 14 mA, 100 ms
Smallest value for current-time product	0.02 mAs
Data according to IEC 60601-1-3	Tube voltage (kV): $\pm 5\%$
Accuracy of loading parameters	Tube current (mA): $\pm 15\%$ (minimum ± 0.1 mA)
	Current-time product (mAs): $\pm 10\%$ (minimum ± 0.2 mAs)
	Pulse duration (ms): $\pm 10\%$ (minimum ± 1 ms)

X-ray tube assembly (Single-tank)

The X-ray tube assembly and image receptor are geometrically aligned such that the emitted radiation cone (radiation axis) is always perpendicular to the plane of the image receptor.

Model	Single tank with single-focus stationary anode tube
Tube material	Tungsten
Nominal focus value	0.6

Tube rated voltage (IEC 613)	110 kV
Single tank rated voltage	110 kV
Anode heat storage capacity (IEC 613)	45,200 J (61,100 HU)
Anode heat dissipation	27,600 J/min (37,300 HU/min)
Anode angle	9°
Inherent filtration tube (IEC 601)	2.8 mm Al with 75 kVp/0.1 mm Cu
Total filtration system (IEC 601)	3.0 mm Al with 75 kVp/0.1 mm Cu
Continuous heat dissipation	80 W
Continuous input power	170 W
Maximum irradiation time	20 min at 600 W / 60 min at 170 W

9.4.3 Image receptor

Flat detector

Cios Fusion utilizes Varian PaxScan detectors. The solid state detector converts the X-rays to visible photons by a scintillator. The visible photons are absorbed by pin photodiodes and converted to electron-hole pairs, which collect on the capacitance of the photodiodes. The pin photodiodes are discharged when the pixel's TFT is turned ON. The charge is collected by an integrating amplifier and converted to a voltage. The signal voltage has a programmable gain applied, depending upon operating conditions. The output voltage is converted to digital data by an analog-to-digital converter.

The following table shows the typical performance test data of the two detector sizes used with the Cios Fusion conducted with the Varian PaxScan solid state detector used with the Cios Fusion.

Detector			Pax Scan 3030XV			Pax Scan 2020XV		
	Unit		Min	Typical	Max	Min	Typical	Max
Technology								
	Scintillator			CsI			CsI	
Geometry								
	Pixel pitch	μm		194			194	
	Active area	mm (V x H)		298 x 298			199 x 199	
	Active area (total) ⁹⁾	Pixels (V x H)		1536 x 1536			1024 x 1024	
	Used binning			1x1, 2x2			1x1 ¹⁰⁾	

9) The useable pixel matrix is reduced on each side by up to 10 pixels relative to the values stated here

10) Although the 2020XV is capable of 2x2 binning as well, this mode is currently not used by Siemens

9 Technical Specifications

Detector			Pax Scan 3030XV			Pax Scan 2020XV		
	Unit		Min	Typical	Max	Min	Typical	Max
Dose characteristics								
Voltage range	kVp		40		150	40		150
Saturation dose (fluoro)	μGy		4.4			8.2		
Saturation dose (single image)	μGy		29			29		
Noise equivalent dose	nGy/ binning			3.5/2x2			3.5/2x2 ¹¹⁾	
Dynamic range ¹²⁾	dB			94	96		94	96
Image quality performance ¹³⁾								
Sensitivity (low-dose fluoro)	LSB/nGy		30	36		17	20	
Electronic noise (low-dose fluoro)	nGy/ binning				0.75/2x2			1.5/1x1
Residual Lag	% in 1st image				5			5
DQE @ 40nGy, RQA5 (30x30 cm², 2x2)								
at 0.05 lp/mm	%			76			n/a	
at 1.0 lp/mm	%			43			n/a	
at 1.29 lp/mm (Nyquist)	%			25			n/a	
DQE @ 80 nGy, RQA5 (20x20 cm², 1x1)								

11) Although the 2020XV is capable of 2x2 binning as well, this mode is currently not used by Siemens

12) Typical dynamic range is defined as $2 \times \log_{10}$ (typical LSB dynamic range) dB, where typical LSB dynamic range = 50000

13) DQE measured according to standard IEC 62220-1

Scattered radiation grid on flat detector (detachable)

C-arm

15) Limited field of view = 20x20 cm² for 3030XV, n/a for 2020XV

9 Technical Specifications

Horizontal movement	20 cm (7.9")
Immersion depth	73 cm (28.7")
Swivel range	± 12°
Vertical movement	43 cm (16.9"), motorized
SID	102 cm (40.2")
Single tank image receptor distance	82 cm (32.3")

Collimator system

Rectangular collimator	For concentric, radiation-free collimation
Slot collimator	For symmetric and asymmetric radiation-free collimation, unlimited rotation

Imaging chain

Flat detector 2020X	1024 (H) x 1024 (V) @ up to 30 fps
Flat detector 3030X	1536 (H) x 1536 (V) @ up to 5 fps
	1024 (H) x 1024 (V) @ up to 30 fps
	768 (H) x 768 (V) @ up to 30 fps
Imaging System	Features, e.g. digital image rotation ± 360°

Monitors

TFT high-brightness color display

Screen diagonal	19" (48 cm)
Image display	1280 x 1024 (pixels)
Brightness (DICOM calibrated)	400 cd/m ²
Maximum brightness, typical	650 cd/m ²
Horizontal/vertical viewing angle	178°/178°
Contrast ratio, typical	900 : 1
Backlight LED technology	

TFT Premium high-brightness color display (optional)

Anti-reflexion coated cover	Glass panel for higher-contrast image quality and easier cleaning
Brightness Uniformity Auto Adjustment	Brightness uniformity is automatically adjusted for a more balanced and uniform image impression
Automatic Brightness Adjustment	Brightness is automatically adjusted based on ambient light

Automatic DICOM Correction	DICOM curves are automatically corrected based on ambient light (This feature greatly speeds up verification of conformity to DIN 6868-157 in affected countries)
Full b/w emulation mode	
Protective collision cover	
Monitor Interface	Input: 1 x DisplayPort
Screen diagonal	19" (48 cm)
Image display	1280 x 1024 (pixels)
Brightness (DICOM calibrated)	500 cd/m ²
Maximum brightness, typical	1000 cd/m ²
Horizontal/vertical viewing angle	170° / 170°
Contrast ratio, typical	1000 : 1
Response time, typical	25 ms
Backlight White LED technology	

Monitor columns

Flex display column (Standard)	Display positioning with rotation angle -30° to 210°; defined lock-in positions at 0° and 180°
Flex Plus display column (optional)	Display positioning with rotation angle -30° to 210°; defined lock-in positions at 0°, 90° and 180°, motorized height adjustment; fold function to protect the displays in park position and during transport

9.4.5 Equipment options

Dose measuring chamber (DAP meter)

Technology	Semiconductor chamber
Active area	Rectangle - 60 mm x 60 mm
DAP resolution	0.01 μGym^2 equivalent to 0.01 cGycm ²
Maximum measurable DAP	3000.00 μGym^2 equivalent to 3000.00 cGycm ²
Reproducibility	< 1% under constant pressure and temperature
Energy range	40 - 150 kV $\pm$ 6%
Attenuation equivalent (Inherent filtration)	< 0.30 mm Al equivalent
Dose rate linearity	Better than $\pm$ 2% for dose area product rate $\geq$ 1 $\mu\text{Gym}^2/\text{s}$, otherwise better than $\pm$ 5%

9 Technical Specifications

Measuring readiness	20 seconds after power on
Work environment	Temperature range: +10 °C to +70 °C Rel. humidity: 15% to 75% (without condensation) Barometric pressure: 700 hPa to 1060 hPa
Precision of the dose measuring chamber (DAP meter)	< 25% (deviation between displayed and measured dose area product)

Automatic dose rate control (ADR)

Precision of automatic dose rate control (ADR)	± 25% (value depends on measurement setup)
--	--

Detector laser light localizer (optional)

Laser class	Class 1 (IEC 60825-1:2007)
Laser type	Semiconductor laser (laser diode)
Wave length	650 nm
Color	Red
Effective power	≤ 3 mW

Single-tank laser light localizer (optional)

Laser class	Class 1 (IEC 60825-1:2007)
Laser type	Semiconductor laser (laser diode)
Wave length	650 nm
Color	Red
Max. power	≤ 3 mW

HD video manager

Outputs	2 x CAT (RJ 45); DVI-D via active CAT/DVI adapter; 1280 x 1024 @ 60 Hz
Inputs	DVI-I (analog or digital) (2x), RGB, YPrPb, Y, S-Video (2x), Composite, SD/HD-SDI (2x), VGA; offers 150 preprogrammable timing settings with up to 1080p (1920 x 1080, 50/60 Hz, interlaced/progressive) including PAL and NTSC

DVI video splitter

Outputs	2 single-link DVI-D; female connectors; up to 1920 x 1200 @ 60 Hz or 1080 p
Conditions for use	The displays or monitors used must be compatible with the SXGA standard 1280 x 1024 at 60 Hz No galvanic separation, therefore the connected monitors must comply with IEC 60601-1-1

SCALANCE W 721 WLAN access point

Data transfer	
Ethernet transfer rate	BASE-T max. 100 MBit/s
Wireless transfer rate	max. 150 MBit/s
Radiated power	typically: 100 mW, max. 200 mW
Supported wireless standards	IEEE 802.11a IEEE 802.11b IEEE 802.11e IEEE 802.11g IEEE 802.11h IEEE 802.11i IEEE 802.11n
Supported standards for power supply	IEEE 802.3at Type 1 (802.3af) and 802.3at Type 2 (power supply via Ethernet)

9.5 List of Standards

Siemens hereby certifies that the Cios Fusion is compliant with the following recognized consensus standards covering electrical and mechanical safety:

Standard organization	Reference number and date	Title of Standard
IEC	60601-1 3rd	Medical Electrical Equipment - Part 1: General Requirements for Safety, 2005
IEC	60601-1-1:2000	Medical electrical equipment - Part 1-1: General requirements for safety - Collateral standard: Safety requirements for medical electrical systems
IEC	60601-1-2 3rd 2007	Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
IEC	60601-1-3: 3rd 2008	Medical electrical equipment - Part 1-3: General requirements for basic safety and essential performance - Collateral Standard: Radiation protection in diagnostic X-ray equipment
IEC	62366: 2007	Medical devices - Application of usability engineering to medical devices

9 Technical Specifications

ISO	14971: 2012	Medical devices - Application of risk management to medical devices
IEC	62304 Ed. 1.0 2006	Medical device software - Software life cycle processes
IEC	60601-2-7 (1998)	Medical electrical equipment - Part 2-7: Particular requirements for the safety of high-voltage generators of diagnostic X-ray generators
IEC	60601-2-28 Edition 2.0 2010-03	Medical electrical equipment - Part 2-28: Particular requirements for the basic safety and essential performance of X-ray tube assemblies for medical diagnosis
IEC	60601-2-32: 1994	Medical electrical equipment - Part 2: Particular requirements for the safety of associated equipment of X-ray equipment
IEC	60601-2-54 2009	Medical electrical equipment - Part 2-54: Particular requirements for the basic safety and essential performance of X-ray equipment for radiography and radioscopy
IEC/PAS	61910-1 2007	IEC PAS 61910-1 First Edition 2007-07, Medical electrical equipment - Radiation dose documentation - Part 1: Equipment for Radiography and Radioscopy
NEMA	PS 3.1 - 3.20 (2011)	Digital Imaging and Communications in Medicine (DICOM) Set
IEC	60825-1:2007	Safety of laser products - Part 1: Equipment classification, and requirements
IEC	IEC 60601-2-43 - Ed. 2.0 2010-03	Medical electrical equipment - Part 2-43: Particular requirements for the safety and essential performance of X-ray equipment for interventional procedures


9.6 Labels

This section provides information on the position of system labels. The numbers and parameters on the labels are examples only. A real system label may contain other numbers and parameters (such as the serial number).


9.6.1 C-arm system

The labels shown below are attached permanently to the following sub-assemblies.

- (1) C-arm travel frame
- (2) Base system (chassis)
- (3) C-arm system control unit
- (4) Flat detector
- (5) Single tank
- (6) Single tank laser and dose measuring chamber

1. Drawer on C-arm travel frame

- (1) Location of labels
- (2) Label behind cover: Loosen knurled screw and pull out drawer

9 Technical Specifications

- (1) Manufacturer's system identification label, approval label

X-ray tube unit

- (1) Supplemental label


- (2) Single tank identification label

- (3) X-ray tube identification label

Flat detector

Second Set of Labels Detector


Collimator

SECOND SET OF LABELS COLLIMATOR

- (1) Supplemental label


(2) Manufacturer's identification label for collimator


(3) Approval-related label


Dose measurement chamber
(optional)

(1) Supplemental label

(2) Identification label

(3) Approval-related label


Single-tank laser light localizer
(optional)

(1) Supplemental label


(2) Identification label

(3) Approval-related label

9 Technical Specifications

General (country-specific)


(1) Identification label (for Brazil only)

2. Base system (chassis)

(1) - (8)

General information

(1) Base system (chassis)

(2) Weight label

(3) Identification label (for Saudi Arabia only)

(4) X-ray control area warning label (for Germany and Austria only)

(5) Warning label (China only)

Hand switch

Model-No 01234567 IP54
Serial-No 10001 Rev01

(6) Hand switch identification label

Standard footswitch

component name
(01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxx MODEL
(422) xxx ORIGIN YY
(21) xxxxx SERIAL IVK

(7) Standard footswitch identification label

Multifunctional footswitch

component name
(01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxx MODEL
(422) xxx ORIGIN YY
(21) xxxxx SERIAL IVK

(8) Multifunctional footswitch identification label

3. C-arm system control unit and remote control unit

(1)

component name
(01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxx MODEL
(422) xxx ORIGIN YY
(21) xxxxx SERIAL IVK

IP65

SIEMENS
RCP-L.L.
(01) 940808907777 GTIN
(240) 1030113 MODEL
(422) 278 ORIGIN DE
(21) 1188 SERIAL IVK
Siemens Healthineers GmbH, Heidenstr. 187
91052 Erlangen, Germany
Made in Germany

IP65

(1) Manufacturer's identification label, protection class label

4. Flat detector, grid, frame, laser light localizer

Flat detector

VARIAN
medical systems
1675 South Pioneer Road
Salt Lake City, UT 84104
U.S.A.
(801) 972-5000

Model: PacScan2020XV
Digital Image Receptor
Part #: 10860783
Serial #: 00000-00
Rev.: AA
MAY 2015

CE
cTUVus
10


component name
(01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxx MODEL
(422) xxx ORIGIN YY
(21) xxxxx SERIAL IVK

(1) Identification labels

THIS RECEPTOR
COMPLIES WITH HHS 21 CFR
SUBCHAPTER J


(3) Approval-related label

Grid



Manufacturer's identification labels

Frame

Manufacturer's identification label, approval-related label

Laser light localizer

Manufacturer's identification label, approval-related label

5. Single tank

Identification label


Warning label (hot surface)

The label is attached to the front, left, and right side of the single tank.

X-ray tube

Identification label


add. filtr. 0.1 Cu

Approval-related label, filter label

Collimator

Identification label



Approval-related label

Generator


Identification label

6. Single tank laser and dose measuring chamber (optional)

(1) Single tank laser

9 Technical Specifications


This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture.
 Manufactured:
 Siemens Healthcare GmbH
 Henkestr. 127, 91052 Erlangen
 Germany

(2) Dose measuring chamber

9.6.2 Monitor trolley

(1) - (9)


This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture.
 Manufactured:
 Siemens Healthcare GmbH
 Henkestr. 127, 91052 Erlangen
 Germany

(1) Manufacturer's identification label, approval-related label

(2) Identification label (for Saudi Arabia only)

(3) Manufacturer's identification label for monitors

WARNING!
 THIS X-RAY UNIT MAY BE DANGEROUS TO PATIENT AND OPERATOR UNLESS SAFE EXPOSURE FACTORS, AND OPERATING INSTRUCTIONS AND MAINTENANCE SCHEDULES ARE OBSERVED.

(4) Warning label

Volt	50 Hz	60 Hz	Long Time	Moment
			Current	Current
100 ~	☐	☐	18 A	29 A
110 ~	☐	☐	17 A	27 A
120 ~	☐	☐	16 A	25 A
127 ~	☐	☐	16 A	24 A
200 ~	☐	☐	10 A	16 A
220 ~	☐	☐	9 A	14 A
230 ~	☐	☐	9 A	14 A
240 ~	☐	☐	9 A	14 A

(5) Supplemental label country-specific connection data

(6) Approval label

(7) Weight label

**This Device is rated for long time operation
Pulsed Radioscopy (duty cycle: 0...99,9%)**

(8) Approval label (not for China)

(9) Warning label read Operator Manual, watch inclines

9 Technical Specifications

10 Options

10.1	System accessories	288
10.1.1	Spacer	289
10.1.2	Grounding cable (optional)	289
10.1.3	Wireless WLAN network connection (optional)	289
	Activating WLAN	290
	Operation via Ethernet cable	290
10.1.4	Multifunctional footswitch (optional)	291
	Selecting the operating mode	291
	Releasing radiation	291
	Storing images (during radiation)	292
	Storing images (after radiation)	292
	Storage for transport	292
10.1.5	Remote control unit (optional)	292
	Attaching the remote control unit	293
	Connecting the remote control unit	294
	Removing the remote control unit and holder	294
	Storing the remote control unit	295
10.1.6	Position memory	296
	Saving a position	296
	Moving to the target position	296
	Deleting the position memory	296
10.1.7	Sterile cover on the C-arm	296
	Applying the sterile cover	297
10.1.8	Sterile cover for the remote control unit	298
10.1.9	Detector laser light localizer (optional)	298
10.1.10	Single-tank laser light localizer (optional)	299
	Protection measures	299
10.1.11	Dose measuring chamber (DAP meter)	299
10.1.12	NaviLink 2D	299
10.1.13	Connection for contrast agent injector	299
	Description of interfaces	300
	Connecting the Metron injector to Cios Fusion	300
	Connecting the Bayer injector to Cios Fusion	301
10.1.14	Video manager	301
10.1.15	DVI video splitter	302
10.1.16	"EMotion" sound system	302
10.1.17	Litho Option	303
10.1.18	LithoReport	303

10.1 System accessories

Anyone who connects additional equipment to the medical device is considered to be configuring the system and is therefore responsible for ensuring that the current system configuration complies with the relevant standards (e.g. system standard IEC/EN 60601-1-1 and/or other applicable standards). In case of questions please speak with your local contact person.

If one of the listed accessory parts requires special operating conditions (e.g. temperature, air pressure, humidity), appropriate attention will be drawn to such in the description. Please follow the Operator Manual provided by the manufacturer.

The following accessories have been approved for use with the Cios Fusion:


Accessories	Part number	Manufacturer
Sterile cover for flat detector	14448832 (FD 20)	Microtec
	14448834 (FD 30)	
Sterile cover for single tank	14448831	Microtec
Sterile cover for remote control unit	14436370	Microtec
Sterile cover with plastic brackets	14448830	Microtec
Grounding cable	14402243	Nicolay Services GmbH
Spacer	14443479	Techno Blech GmbH
Wireless network connection WLAN	14436364	Siemens Healthcare GmbH
Multifunctional foot-switch	14436357	Steute
Remote control unit and holder	14443478	Penta
Additional holder for remote control unit	14436360	Siemens Healthcare GmbH
Integrated laser light localizer for flat detector	14436342	Siemens Healthcare GmbH
Single tank laser light localizer	14436343	Siemens Healthcare GmbH
Non-sterile cover for footswitch	14436371	Microtec

Accessories	Part number	Manufacturer
Non-sterile cover for multifunctional foot-switch	14436372	Microtec
Arterion Mark7 Interface cable Cios	11154936	Siemens Healthcare GmbH
Mark V ProVis Interface cable Cios	11154943	Siemens Healthcare GmbH

10.1.1 Spacer

The distance between the source and tube assembly cover (shortest possible source-skin distance) is ≥ 200 mm with the standard system (according to IEC 601-1-3).

Country-specific regulations may require a larger source-skin distance (≥ 300 mm according to DHHS 21 CFR). This is achieved by attaching an additional spacer to the C-arm system.

- ◆ If this source-skin distance is too large for special examinations, the spacer can be removed by taking off the screws.


The spacer must be reattached after this type of examination to ensure the reduction in skin dose resulting from a greater source-skin distance.

10.1.2 Grounding cable (optional)

An optional grounding cable is available for equipotential bonding in accordance with DIN 57107/VDE 0107 for rooms of Application Group 2E (cardiac catheterization).

10.1.3 Wireless WLAN network connection (optional)

The WLAN option enables a wireless connection between the Cios Fusion and a (hospital) network. The dedicated WLAN unit is integrated into the monitor trolley housing and is completely set up. Its On/Off switch is found on the side of the monitor trolley.

(1) WLAN On/Off switch

Activating WLAN

- ◆ To activate WLAN, use the switch on the side of the monitor trolley.
The WLAN switch illuminates.

Once the WLAN unit is switched on, it takes less than a minute for the WLAN connection to be functional.

CAUTION

Interference.

Effect on persons with implants.

- ◆ Maintain a minimum distance in accordance with IEC 60601-1-2. If this is not possible, switch off the WLAN client.

The coexistence of wireless sources operated in the same environment may affect the quality of data transfer.

The transmission speed is reduced if additional WLAN clients are operated simultaneously in the same environment.

For reasons of data security, it is recommended to encrypt the data in the WLAN according to WPA or WPA2.

In case the wireless connection is interrupted, impaired or overloaded, it is recommended to keep an Ethernet cable on hand, provided there are connection options nearby.

If WLAN is activated, no devices sensitive to high frequencies should be placed on top of the unit. Otherwise, damage to such devices or impairment to WLAN functionality could result.

Operation via Ethernet cable

When operating the Cios Fusion via Ethernet cable, the WLAN at the monitor trolley needs to be deactivated prior to switching on the system.

Simultaneous operation via Ethernet cable and WLAN is not possible and should be avoided by all means.

After inadvertent simultaneous operation or after a change between both options, the system must be switched off and on again.

10.1.4 Multifunctional footswitch (optional)

Selecting the operating mode

The MODE (operating mode) button of the multifunctional footswitch allows you to select one of the operating modes: Single image, fluoro, SUB, Roadmap or a video source (depending on the configuration level).

● MODE

- ◆ If necessary, press this button on the multifunctional footswitch several times.

The selected operating mode is indicated in the **Examination** task card and on the control panel of the C-arm system.

The set video source alters the image content of the right-hand monitor on the monitor trolley.

A confirmation sound acknowledges successful operating mode/video source switchover. You can then release radiation again right away.

Releasing radiation

The left pedal is always used to activate fluoroscopy (CFC) (Preferred setting).

The right pedal is used to activate the currently selected operating mode.

Exception: If Fluoroscopy (CFC) mode is selected, the right pedal is assigned the single image operating mode.

The functionality of the pedals can optionally be changed.

- ◆ Keep the foot pedal pressed during radiation release.

Storing images (during radiation)

- ◆ Press the save button on the multifunctional footswitch during radiation.
The image currently generated and displayed is saved in the local database.

Storing images (after radiation)

- ◆ Press the save button on the multifunctional footswitch.
Holding the key for < 2 seconds: saves the image last recorded (LIH).
Holding the key for > 2 seconds: saves the scene last recorded (LSH).


The images are stored in the local database. They are transferred from the left monitor to the right monitor unless the **Hold reference** function is active.

Storage for transport

To transport the multifunctional footswitch properly, proceed as follows:

- ◆ Insert the multifunctional footswitch and wrapped cable into the holder provided on the front of the C-arm system.

10.1.5 Remote control unit (optional)


The remote control unit allows the C-arm system to be operated in the OR directly from the patient table. It is used in the same way as the control unit on the C-arm system.

CAUTION

The surfaces of the remote control unit become hot.

Risk of burns!

- ◆ Make sure that the patient does not touch the remote control unit for more than 10 minutes or that after 10 minutes can no longer come in contact with it.
- ◆ As a rule, contact between the patient and the remote control unit should be avoided.

(1) Integrated EMERGENCY STOP button

(2) Control panel

Attaching the remote control unit

CAUTION


Remote control unit falls down.

Risk of injury/crushing!


- ◆ It is recommended to keep the remote control unit in the holder provided. If a holder is not available, place the remote control unit on the instrument table.

Attach the remote control unit to the patient table railing as follows:


- 1 Attach the holder to the table railing (1) and tighten the holder screw (2).

- 2 Place the patient sterile cover over the holder.

- 3 If necessary, pack the remote control unit in the sterile cover and place it on the magnetic holding plate.


CAUTION

Remote control unit covered by unsuitable covers.

Overheating of the remote control unit!

◆ Only use the sterile covers specified!

Connecting the remote control unit

(1) Round socket for connecting the remote control unit

- ◆ Plug the connector of the remote control unit into the round socket on the C-arm system, which is located next to the central plug.


The remote control unit is ready for operation and can be used in parallel with the control panel on the C-arm system.


The remote control unit can also be connected to other C-arm systems. In this case a message is displayed that you have to confirm.

Removing the remote control unit and holder

- 1 Unplug the connector of the remote control unit from the round socket on the C-arm system.
- 2 Remove the remote control unit from the magnet by slightly tilting it to the side.

- 3 Remove the holder by first turning the screw (1) and then removing the holder upward (2).

Storing the remote control unit

CAUTION



Improper positioning of the remote control unit.

Risk of contamination.

- ◆ Do not place the remote control unit on the patient.

To store the remote control unit, proceed as follows:

- 1 Remove the remote control unit together with the holding plate from the table railing and wrap the cable around the back.

- 2 Attach the remote control unit together with the holding plate to the monitor trolley railing (see figures 1 and 2).

10.1.6 Position memory

With this option, the orbital and angular positions of the C-arm will be displayed on the control panels and the left monitor. You can save an orbital and an angular position in the position memory. After intermediate changes of acquisition position, this enables you to move the C-arm back to the original position.

Saving a position

- 1 Perform angular and orbital C-arm movement to a target position you want to save.

- 2 Press this button on the control unit.

The current orbital and angular positions are saved.


Moving to the target position

- ◆ Move the C-Arm toward the target position.

The required direction of movement and progress toward reaching the target is indicated by corresponding icons on the control panel, see the following example illustration:

Icon	Meaning
	Deviation from target position > 5°
	Deviation from target position < 5°
	Deviation from target position < 1°
	Target position reached precisely

On the left monitor, an arrow of the same color shows the deviation and the required direction of movement during the examination.

Deleting the position memory


- ◆ Press the button in the activated state.

The position memory is deleted.

10.1.7 Sterile cover on the C-arm

To protect it against contamination, the C-arm including the flat detector and X-ray tube assembly are completely covered with a three-piece sterile disposable cover of transparent polythene sheet.

C-arm completely covered

CAUTION

Improper attachment of the sterile cover on the C-arm.

Risk of infection!

♦ Follow the instructions for handling the sterile covers.


Always use sterile covers during interventions to protect the flat panel detector against the ingress of liquids!


Applying the sterile cover

Follow these instructions to wrap the C-arm in the sterile cover:

- 1 Attach the disposable sheet to the C-arm using plastic clips.


- 2 Pull the plastic cover over the single tank.


The plastic cover is fixed in place with an elastic cord.

- 3 Pull the other plastic cover over the flat detector.


- 4 Secure the plastic cover over the handle in the locations provided.

The C-arm is now covered completely.

10.1.8 Sterile cover for the remote control unit

- ◆ Insert the remote control unit in the tube-like plastic sack and close it carefully.

10.1.9 Detector laser light localizer (optional)

To project the target crosshairs, two Class 1 lasers are used whose exit windows are arranged in the detector cover.

10.1.10 Single-tank laser light localizer (optional)

To project the target crosshairs, two Class 1 lasers are used whose exit windows are arranged in the single tank cover.

Protection measures


Brief exposure to a class 1 laser is safe for the eye.

The aversion response and the blinking reflex will usually protect the eye.

In Germany, the operator is responsible for ensuring that the user has been instructed in the use of the laser. Outside Germany, the relevant laws and regulations regarding the use of Class 1 lasers must be complied with.

10.1.11 Dose measuring chamber (DAP meter)

If available, this option is integrated in the single tank of the Cios Fusion. The DAP meter determines the dose area product (DAP) and air kerma values (→ Page 238 *Dose and consistency test*). The cumulative patient data are displayed on the monitor trolley during the examination (→ Page 126 *Radiation information*).

Dose measuring chamber (DAP meter), integrated in the single tank

10.1.12 NaviLink 2D

If the **NaviLink 2D** option is available in the system and communication between the navigation system and the C-arm is activated, you can see a **NaviLink 2D** icon in the status bar (bottom left) of the left monitor.

The navigation system sends a message to the C-arm system signaling readiness for receiving images. The C-arm system receives the message, identifies the sending navigation system and responds. The communication connection between the two systems is activated and the **NaviLink 2D** icon displays. As soon as the communication connection is terminated, the **NaviLink 2D** icon disappears again.

10.1.13 Connection for contrast agent injector

A connection (unidirectional trigger output) for a contrast agent injector is possible as an option.

The use of a contrast agent injector can result in an improved image quality during subtraction angiography, since the stable flow rate of the contrast agent ensures homogeneous contrast enhancement during the entire examination. This can also reduce the patient's exposure to contrast agent and radiation.

Refer to the operator manual from the respective manufacturer for the operation and start-up of the released contrast agent injector.

Description of interfaces

If the injector connected to this interface is in the "armed" state, Cios Fusion provides the radiation on signal. This signal can be used by the injector to start the contrast agent injection.

Since the optimal point in time for the injection of contrast agent is right after the subtraction mask has been created, the delay on the injector shall be set equal to the mask creation time. The mask creation time is directly proportional to the currently set pulse rate. The mask creation time in seconds is displayed on the left monitor close to the subtraction mask progression bar. It is recommended to set the delay on the injector equal to the mask creation time.

CAUTION

Radiation release without taking the contrast agent injector status into account.

Risk of unnecessary radiation exposure.


- ◆ Therefore, make sure that the injector is ready to administer contrast agent before releasing radiation.

There are two injectors available for the Cios Fusion:

- Metron injector
- Bayer injector


Connecting the Metron injector to Cios Fusion

- 1 Attach the contrast agent injector interface electronics housing to the Cios Fusion electronics unit using the attachment kit provided.

- 2 Insert the interface electronic housing connector into the corresponding Cios Fusion electronics unit socket.


Connecting the Bayer injector to Cios Fusion

A system-specific connection cable is used for wired systems to connect the injector.

- ◆ Connect the Bayer injector interface cable to the corresponding Cios Fusion electronics unit socket.


10.1.14 Video manager

Cios Fusion has an additional HD video manager option for emitting video signals, e.g. on external monitors, or for feeding in video signals for display on the reference monitor (monitor B), e.g. for endoscopy images.


Only medical grade monitors must be connected.

10 Options

The video signal feed interfaces are located on the back of the monitor trolley (see left graphic).

To establish the connection to an external monitor, proceed as follows.

- 1 Insert and secure the receiver module (1) in the DVI interface of your external monitor.
- 2 Using a CAT cable (2), connect the receiver to the respective socket A (live monitor) or socket B (reference monitor) on the monitor trolley.
- 3 Insert the 5 V connector of the power supply (3) into the receiver.
- 4 Plug the corresponding power plug into the power outlet.

The connection to the external monitor is now operational.

10.1.15 DVI video splitter

Cios Fusion has an additional DVI video splitter option available for outputting image signals to an external live monitor (monitor A) and an external reference monitor (monitor B) via the DVI connection.

Only medical grade monitors must be connected.

- ◆ To connect an external monitor, insert the DVI connection cable into the respective interface for monitor A or B on the monitor trolley.

10.1.16 "EMotion" sound system

The EMotion integrated sound system is available as an option for Cios Fusion, with an interface for operating external (MP3) audio equipment via Aux-in jack (Aux input) including amplifiers with digital sound processor (DSP) and 2-way stereo loudspeaker systems for high quality sound.


Do not connect a device with its own power cord to the Aux input, use battery-operated devices only.


10.1.17 Litho Option

The Cios Fusion "Litho" allows the C-arm system to be combined with a MODULARIS Variostar ESWL unit. This is accomplished via two different connector ports.

The MODULARIS Variostar hand control unit can be connected to the port on the front of the C-arm system and used to create and save X-ray images:

The hand control unit of the MODULARIS Variostar can be connected to the second port to control the motorized angulation of the C-arm:

Instructions for the operation and initial startup of the Cios Fusion Litho with MODULARIS Variostar can be found in the appropriate Operator Manual.

10.1.18 LithoReport

The LithoReport option provides the Cios Fusion Litho with a comprehensive tool for the professional management of patient data, including diagnostic, therapeutic, shock wave, and continued care data.

Instructions for the operation and initial startup of LithoReport can be found in the appropriate Operator Manual.

This page has been intentionally left blank.

This page has been intentionally left blank.

10 Options

This page has been intentionally left blank.

Manufacturer's notes:

This product is provided with a CE marking in accordance with the regulations stated in Annex II of Directive 93/42/EEC of June 14th, 1993 concerning medical devices, as well as directive 2011/65/EU of June 8, 2011 regarding the restricted use of certain hazardous substances in electrical and electronic equipment.

The CE marking applies only to medical devices which have been placed on the market according to the above-mentioned EC Directive.

Unauthorized changes to this product invalidate this declaration.

Original language: English

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Legal Manufacturer

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

Print No. XPR2-370G.620.01.02.02 | © Siemens Healthcare GmbH, 2017


www.siemens.com/healthcare

Deed No. N 0147 /2022 Fr

I, Marcel Neumair, notary in Forchheim, notarize that the document overleaf was signed in my presence by

Mr. Ali-Ihsan **Can**,
born on December 11th, 1963,
with business address Siemens Healthcare GmbH,
Siemensstr. 1, 91301 Forchheim, Germany,
known to me, the notary, by previous notarization.

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The screenshot shows a window titled "Query for retrieve". It is divided into two main sections: "PATIENT" and "STUDY".

PATIENT section:

- Last Name: [Text input field]
- First Name: [Text input field]
- Patient ID: [Text input field]
- Date of Birth: [D:] [M:] [Y:]
- Sex: ☐ Male ☐ Female ☐ Other

STUDY section:

- Accession No: [Text input field]
- Study ID: [Text input field]
- Node: [Dropdown menu showing "20"]

Results section:

Below the input fields is a table with the following columns: Patient Name, Patient ID, Date of Birth, Sex, Study Description, Date and, Access, and Ins. The table contains several rows of data.

At the bottom of the window is a button labeled "Studies".

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Name	Karla Flesner-Oeste
Department	SHS AT QT RA
Mobile	+49 (1525) 9518360
E-mail	karla.flesner-oeste@siemens-healthineers.com
Date	January 24, 2022


Declaration Letter Operator Manual Addendum

We, Siemens Healthcare GmbH, hereby confirm that the following Operator Manual Addendum for
Cios Fusion has been printed by the manufacturer in Germany:

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www.siemens.com/healthineers

Cios Fusion Important additional information

Operator Manual

This device bears a CE mark in accordance with the provisions of Council Directive 93/42/EEC of June 14, 1993 concerning medical devices and the Council Directive 2011/65/EU of June 08, 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

The CE marking applies only to Medical Devices which have been put on the market according to the above-mentioned EC Directives.

Unauthorized changes to this product are not covered by the CE mark and the related Declaration of Conformity.

Original language: English

Legal Manufacturer

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
[siemens.com/healthcare](https://www.siemens.com/healthcare)

www.siemens.com/healthineers

Operator Manual

Important additional information

Use of the product	3
Contraindications	3
Undesireable side effects	3
Residual risk	3
Safety	5
Security	5
Malware protection	5
System Description	5
Reconnecting the monitor trolley	5
Radiation release	5
Image display	5
Function displays	5
C-arm system park position	6
Examination	6
Selecting the application	6
Temperature monitoring	7
Radiation information	7
Post-processing	8
Reviewing imported scenes	8
Technical Specification	9
Only for Ukraine	9
Only for Australia	9
Precautions for EMC – 3rd Edition	10
EMC Regulations - 4th Edition	10
Clarification on requirement compliance with regard to interference emission and immunity	10
Clarification regarding EMC environment	12
Operation directly next to or stacked with other devices	12
Listed cables, transformers, and accessories	12
Use of other accessories, cables, and transformers	12
Portable RF devices	12
Note regarding Class A devices according to CISPR 11	12
Reduction or loss of critical features due to electromagnetic interference	13
Instructions for maintaining the basic safety and critical features in reference to electromagnetic interference for the expected useful life of the system	13

Important additional information

Table of Contents

Use of the device near transmitters	13
List of standards	14
Image quality performance	14
Dosimetric Information: valid until VA20G delivery state	14
Radioscopy, Detector size 20x20	14
Radioscopy, Detector size 30x30	16
Radiography, Detector size 20x20	17
Radiography, Detector size 30x30	18
Dosimetric Information: valid from VA20H delivery state	20
Radioscopy, Detector size 20x20	20
Radioscopy, Detector size 30x30	21
.....	21
Radiography, Detector size 20x20	23
Radiography, Detector size 30x30	24
Options	26
Video manager	26
List of Warnings and Cautions	27

Important additional information

Please note, that the following descriptions, notes, and technical data replace the respective text written in the system Operator Manual Cios Fusion.

Use of the product

Contraindications

There are no known contraindications for this product at present. The physician makes the final decision regarding the use of the medical device for a specific application based on his or her medical knowledge, while weighing the inherent risks for the respective cases.

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In combination with the appropriate care in the clinical workflow and informed decisions considering also the clinical benefit of radiation, the above mechanisms should prevent deterministic dose effects to the patient (erythema, cataracts, permanent epilation and delayed skin necrosis). Nevertheless, such injuries cannot be completely excluded. In rare cases, the dose management mechanisms may be not be operating correctly. These issues are detectable for the user especially within a Quality Assurance framework. Some clinical applications using the Cios Fusion involve injecting contrast medium into the patient's cardiovascular system or other body orifices. Functional disturbances in the Cios Fusion may affect the imaging workflow resulting in the additional administration of contrast medium. This may lead to risks and side-effects in accordance with the information provided by the manufacturer.

Safety

Security

The preinstalled security package offers protection against cyberattacks, viruses, malware, and other malicious software. It ensures that only trustworthy applications can be run on the system. It blocks unauthorized access attempts, protects against network threats and infected USB sticks, and offers control over who uses which applications and when.

Malware protection

Malware protection is based on a whitelisting solution (McAfee Embedded Control) which is pre-installed and pre-configured on the system. This ensures that only manufacturer-authorized software or updates can be installed and executed.

System Description

Reconnecting the monitor trolley

Before reconnecting the monitor trolley to the base unit, make sure that the image system has fully booted.

Connecting the monitor trolley before or just as the starting screen displays can cause the system to freeze. In this case, the system has to be restarted.

When reconnecting, always wait until the connection is fully established to ensure full system functionality.

Radiation release

When pressing the hand switch or foot switch, an audible signal is heard (if configured) until radiation is applied. The length of the tone can vary depending on the operating mode.

Image display

Function displays

Important additional information

The displays for the DAP meter can be specifically configured by service, for example:

"Gycm2" "Gy*cm²"

"dGycm2" "dGy*cm²"

"cGycm2" "cGy*cm²"


"microGym2" "μGy*m²" (DEFAULT)

"mGycm2" "mGy*cm²"

C-arm system park position

When parking the C-arm system on an inclined surface, proceed as follows:

- ◆ Lock the brakes on the C-arm system.
 - The red indicator appears (brakes locked).

- ◆ Set the steering lever to lock position for transverse travel to the right or left, so that the wheels are braked in transverse position.

Examination

Selecting the application

In the PEX Editor, if in one of the applications a curve with the extension HL (High Level) was programmed for a "Low" or "Medium" dose level, or if an ESU (electric noise suppression) characteristic was programmed and the system is not equipped with ESU, an error message will be displayed when this function is selected. In this case, however, radiation can still be released.

Temperature monitoring

Where necessary, fluoroscopy power is limited to prevent overheating.

Storage capacity

The storage capacity refers only to image storage, not to the overall storage capacity of the hard drive.

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00:00:11
hh:mm:ss

The cumulative dose area product for the current patient is displayed if the dose measurement chamber is installed.

55
μGym²

As an alternative: Display of the air kerma value (2) and cumulative air kerma value (1). These data appear on the left monitor.

(1) — 0.4 mGy
00:00:29 hh:mm:ss
(2) — 0.8 mGy/min

The reference location for determining the air kerma values with the dose measurement chamber is 30 cm in front of the detector input. The reason for this convention is that in typical applications the object to be examined is located approximately 30 cm in front of the flat detector.

Depending on country-specific regulations, you can have Siemens Customer Service change the display so that it indicates the air kerma and cumulative air kerma

value instead of the dose area product.

Post-processing


Reviewing imported scenes

If you are reviewing an imported scene with a high frame rate that originated from another modality (such as a CT scene), you may not be able to stop the playback, in which case the system has to be shut down to end the infinite loop.

Always review scenes imported from other modalities at a normal frame rate.


Technical Specification

The detector manufacturer name has changed from VARIAN to VAREX.


The corresponding labels of the flat detector

Only for Ukraine

Only for Australia

Precautions for EMC – 3rd Edition

According to the IEC 60601-2-54 Standard the product provides following Essential Performance

- ❑ Accuracy of Loading Factors
- ❑ Reproducibility of the Radiation output
- ❑ Automatic Control System
- ❑ Imaging performance

According to the IEC 60601-2-43 Standard the product provides following Essential Performance

- ❑ Recovery management
- ❑ Radiation dose documentation

EMC Regulations - 4th Edition

Clarification on requirement compliance with regard to interference emission and immunity

Cios Fusion meets the following test level of standard IEC 60601-1-2:2014

Interference emission tests:	EMC base standard or test procedure:	Compliance level:
Conducted or radiated RF interference	CISPR 11	Group 1, Class A

Immunity testing:	EMC base standard or test procedure:	Compliance level:
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2, 4, 8, 15 kV air discharge
Rapid electric transients/bursts	IEC 61000-4-4	± 2 kV for power supply lines; 100 kHz repetition rate ± 1 kV for input and output cables; 100 kHz repetition rate
Surges	IEC 61000-4-5	± 1; 2 kV conductor-ground ± 1 kV conductor-conductor

Immunity testing:	EMC base standard or test procedure:	Compliance level:
Voltage dips on power lines	IEC 61000-4-11	0% UT for 0.5 cycles (1 phase) 0% UT for 1 cycle 70% UT for 25/30 cycles (50/60 Hz)
Voltage interrupts on power lines	IEC 61000-4-11	0% UT; 250/300 cycles
Power frequency (50 or 60 Hz) magnetic field	IEC 61000-4-8	30 A/m
Conducted interference induced by RF fields	IEC 61000-4-6	150 kHz to 80 MHz: 3 Vrms 6 Vrms in ISM bands between 0.15 MHz and 80 MHz
Radiated electromagnetic RF fields	IEC 61000-4-3	80 MHz – 2.7 GHz: 3 V/m
Induction fields of wireless RF communication devices	IEC 61000-4-3	385 MHz: 27 V/m; pulse modulation 18 Hz 450 MHz: 28 V/m; frequency modulation ± 5 kHz deviation 1 kHz sine 710, 745, 780 MHz: 9 V/m; pulse modulation 217 Hz 810, 870, 930 MHz: 28 V/m; pulse modulation 18 Hz 1720, 1845, 1970 MHz: 28 V/m; pulse modulation 217 Hz 2450 MHz: 28 V/m; pulse modulation 217 Hz 5240, 5500, 5785 MHz: 9 V/m; pulse modulation 217 Hz

Clarification regarding EMC environment

The Cios Fusion is suitable for professional medical care facilities in accordance with IEC 60601-1-2, 4th edition.

Operation directly next to or stacked with other devices

Avoid using this device directly next to or stacked with other devices, as this may lead to improper functioning of the device. If such use is unavoidable, observe whether the affected devices function normally.

Listed cables, transformers, and accessories

Fixed product cables that cannot be removed by the user, are not listed. These cables are part of the product and were taken into account during the EMC measurements. The product is not fully functional without these cables. Deviations or supplements to this document can be found in the product-specific documentation.

Use of other accessories, cables, and transformers

The use of accessories, transformers, and cables other than those specified or delivered by the manufacturer of the device can lead to increased electromagnetic interference or reduced electromagnetic immunity of this device, with the result that the device does not properly function.

Portable RF devices

Portable RF communication devices (including peripheral devices such as antenna cables and external antennas) must maintain a minimum distance of 30 cm (12 inches) from the system (including the cable specified by the manufacturer). Otherwise, reduced device performance is possible.

Note regarding Class A devices according to CISPR 11

Due to its emission characteristics, this device is suitable for operation in industrial areas and hospitals (Class A according to CISPR 11). When operated in a residential area (where generally a CISPR 11 classification as a Class B device is required), the system may not provide adequate protection against RF communication services. Where needed, the user will have to take measures to reduce interference, such as a modified layout or a change in location for the system.

Reduction or loss of critical features due to electromagnetic interference

Caution

Electromagnetic system interference.

Reduction in image quality or the accuracy of radiation parameters.

- ◆ Remove the source of interference (such as a cell phone) from the area.
- ◆ Shut down the system and switch it back on.
- ◆ Watch for error messages during the restart.

Instructions for maintaining the basic safety and critical features in reference to electromagnetic interference for the expected useful life of the system

Perform an annual visual inspection of the accessible plug connectors and screws to ensure there has been no degradation in electrical connections. If there has been degradation, contact Siemens Customer service to repair or replace the corresponding part. This ensures that the design measures regarding electromagnetic compatibility remain effective for the entire useful life of the device.

Use of the device near transmitters

This device has been checked for immunity against RF interference at select frequencies only. The use of transmitters in the immediate vicinity that utilize other frequencies could result in the device not functioning as intended.

List of standards

The standards list 60601-2-7 (1998) is obsolete and has been fully integrated into the standards list 60601-2-54 (2015). It is therefore no longer part of the operator manual.

Image quality performance

Detector		PaxScan 3030XV	PaxScan 2020XV
	Unit	Typical	Typical
MTF (1x1 binning)			
at 1.0 lp/mm	%	55	55
at 2.0 lp/mm	%	25	25
at 2.6 lp/mm (Nyquist)	%	16	16

Dosimetric Information: valid until VA20G delivery state

Radioscopy, Detector size 20x20

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate at patient entrance reference point [mGy/min]
Intervention, Cardiac Interv Zoom 0, 20 x 20	low	80	11.7	7.1	15	3.8
	medium	80	2.8	cont.	cont.	9.2
	high	80	4.1	cont.	cont.	13.4
Intervention, Cardiac Interv Zoom 1, 15 x 15	low	80	14.7	7.1	15	4.9
	medium	80	3.8	cont.	cont.	11.7
	high	80	5.5	cont.	cont.	17.2

Important additional information

Intervention, Cardiac Interv Zoom 2, 10 x 10	low	80	20.9	7.1	15	7.1
	medium	80	5.8	cont.	cont.	19.3
	high	80	8.3	cont.	cont.	26.2
Intervention, Gastro Interv Zoom 0, 20 x 20	low	82	10.0	7.1	15	3.5
	medium	91	1.5	cont.	cont.	6.6
	high	94	2.0	cont.	cont.	9.0
Intervention, Vasc Interv Zoom 0, 20 x 20	low	63	24.0	12.4	7.5	3.8
	medium	65	23.6	13.1	15	9.0
	high	66	23.5	17.5	15	12.4
Cardiac, Card Standard Zoom 0, 20 x 20	low	80	2.0	cont.	cont.	6.3
	medium	80	2.9	cont.	cont.	9.1
	high	80	4.2	cont.	cont.	13.5
Gastro, Gastro Standard Zoom 0, 20 x 20	low	88	1.3	cont.	cont.	5.3
	medium	91	1.5	cont.	cont.	6.6
	high	94	2.0	cont.	cont.	9.0
Vascular, Vasc Standard Zoom 0, 20 x 20	low	63	23.8	8.9	15	5.4
	medium	65	23.7	11.7	15	7.7
	high	67	23.5	15.8	15	11.5
Urology, Uro Standard Zoom 0, 20 x 20	low	64	23.8	10.3	15	6.5
	medium	65	23.7	13.1	15	8.9
	high	67	23.5	16.3	15	11.9
Ortho/Trauma, Ortho Standard Zoom 0, 20 x 20	low	64	23.8	10.3	15	6.5
	medium	65	23.6	13.1	15	9.0
	high	67	23.5	16.3	15	11.9
General, Standard Zoom 0, 20 x 20	low	64	23.8	10.3	15	6.5
	medium	65	23.6	13.1	15	8.9
	high	67	23.5	16.3	15	11.9
General, Standard HL Zoom 0, 20 x 20	low	80	10.8	7.1	15	3.6
	medium	80	18.5	7.1	15	5.8
	high	80	23.6	13.0	15	14.6
Intervention, Cardiac Interv (high dose) Zoom 0, 20 x 20	high dose	110	8.8	cont.	cont.	54.5
General, Pediatrics without scattered radiation grid Zoom 0, 20 x 20	low	71	6.3	7.1	15	1.5
	medium	74	6.6	7.1	15	1.8
	high	78	7.0	7.1	15	2.2

Important additional information

Radioscopy, Detector size 30x30

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate at patient entrance reference point [mGy/min]
Intervention, Cardiac Interv Zoom 0, 30 x 30	low	79	8.4	7.1	15	2.2
	medium	80	1.6	cont.	cont.	5.0
	high	80	2.3	cont.	cont.	7.5
Intervention, Cardiac Interv Zoom 1, 20 x 20	low	80	11.3	7.1	15	3.7
	medium	80	2.8	cont.	cont.	9.2
	high	80	4.1	cont.	cont.	13.5
Intervention, Cardiac Interv Zoom 2, 15 x 15	low	80	14.1	7.1	15	4.9
	medium	80	3.8	cont.	cont.	11.9
	high	80	5.5	cont.	cont.	17.8
Intervention, Gastro Interv Zoom 0, 30 x 30	low	78	6.8	7.1	15	2.2
	medium	87	1.0	cont.	cont.	4.0
	high	90	1.4	cont.	cont.	5.5
Intervention, Vasc Interv Zoom 0, 30 x 30	low	61	24.0	8.6	7.5	2.3
	medium	63	24.0	9.2	15	5.5
	high	64	23.5	12.4	15	7.6
Cardiac, Card Standard Zoom 0, 30 x 30	low	80	1.1	cont.	cont.	3.4
	medium	80	1.6	cont.	cont.	5.0
	high	80	2.4	cont.	cont.	7.5
Gastro, Gastro Standard Zoom 0, 30 x 30	low	85	0.8	cont.	cont.	2.9
	medium	87	1.0	cont.	cont.	4.0
	high	90	1.4	cont.	cont.	5.5
Vascular, Vasc Standard Zoom 0, 30 x 30	low	62	23.7	7.6	15	4.2
	medium	63	23.6	9.2	15	5.5
	high	64	23.5	11.5	15	7.2
Urology, Uro Standard Zoom 0, 30 x 30	low	62	23.7	7.6	15	4.2
	medium	63	23.6	9.0	15	5.4
	high	64	23.5	11.5	15	7.2
Ortho/Trauma, Ortho Standard Zoom 0, 30 x 30	low	62	23.7	7.6	15	4.2
	medium	63	23.6	9.2	15	5.5
	high	64	23.5	11.5	15	7.2

General, Standard Zoom 0, 30 x 30	low	62	23.7	7.6	15	4.2
	medium	63	23.6	9.0	15	5.3
	high	64	23.5	11.5	15	7.3
General, Standard HL Zoom 0, 30 x 30	low	78	6.8	7.1	15	2.2
	medium	80	10.9	7.1	15	3.2
	high	80	23.6	7.7	15	8.1
General, Standard HL (high dose) Zoom 0, 30 x 30	high dose	110	8.8	cont.	cont.	50.8
General, Pediatrics without scattered radiation grid Zoom 0, 30 x 30	low	65	5.6	7.1	15	1.1
	medium	68	5.9	7.1	15	1.3
	high	72	6.2	7.1	15	1.6

Radiography, Detector size 20x20

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Charge [mAs]	Dose at patient entrance reference point [mGy]
Intervention, Cardiac Interv Zoom 0, 20 x 20	low	80	1.34	0.1
	medium	80	4.07	0.3
	high	80	7.91	0.5
Intervention, Cardiac Interv Zoom 1, 15 x 15	low	80	1.76	0.1
	medium	80	5.32	0.3
	high	80	10.44	0.6
Intervention, Cardiac Interv Zoom 2, 10 x 10	low	80	2.64	0.2
	medium	80	8.14	0.4
	high	80	15.99	0.8
Intervention, Gastro Interv Zoom 0, 20 x 20	low	92	0.68	< 0.1
	medium	95	1.87	0.2
	high	94	3.79	0.3
Intervention, Vasc Interv Zoom 0, 20 x 20	low	74	1.95	0.1
	medium	76	5.24	0.3
	high	76	10.44	0.5
Cardiac, Card Standard Zoom 0, 20 x 20	low	80	1.30	0.1
	medium	96	1.83	0.2
	high	110	1.57	0.2

Important additional information

Gastro, Gastro Standard Zoom 0, 20 x 20	low	92	0.65	<0.1
	medium	98	1.62	0.2
	high	110	1.58	0.2
Vascular, Vasc Standard Zoom 0, 20 x 20	low	73	2.12	0.1
	medium	94	1.88	0.2
	high	110	1.57	0.2
Urology, Uro Standard Zoom 0, 20 x 20	low	73	2.11	0.1
	medium	94	1.88	0.2
	high	110	1.57	0.2
Ortho/Trauma, Ortho Standard Zoom 0, 20 x 20	low	73	2.11	0.1
	medium	94	1.88	0.2
	high	110	1.57	0.2
General, Standard Zoom 0, 20 x 20	low	73	2.11	0.1
	medium	94	1.88	0.2
	high	110	1.57	0.2
General, Standard HL Zoom 0, 20 x 20	low	80	1.30	0.1
	medium	96	1.84	0.2
	high	110	1.57	0.2
General, Standard HL (high dose) Zoom 0, 20 x 20	high dose	110	14.08	1.4
General, Pediatrics without scattered radiation grid Zoom 0, 20 x 20	low	87	0.24	<0.1
	medium	92	0.59	<0.1
	high	95	1.06	0.1

Radiography, Detector size 30x30

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Charge [mAs]	Dose at patient entrance reference point [mGy]
Intervention, Cardiac Interv Zoom 0, 30 x 30	low	80	0.73	< 0.1
	medium	80	2.46	0.1
	high	80	4.54	0.2
Intervention, Cardiac Interv Zoom 1, 20 x 20	low	80	1.30	0.1
	medium	80	3.92	0.2
	high	80	7.52	0.4

Important additional information

Intervention, Cardiac Interv Zoom 2, 15 x 15	low	80	1.73	0.1
	medium	80	5.22	0.3
	high	80	10.35	0.5
Intervention, Gastro Interv Zoom 0, 30 x 30	low	88	0.47	< 0.1
	medium	91	1.38	0.1
	high	91	2.66	0.2
Intervention, Vasc Interv Zoom 0, 30 x 30	low	72	1.29	<0.1
	medium	73	3.38	0.2
	high	73	7.21	0.3
Cardiac, Card Standard Zoom 0, 30 x 30	low	80	0.75	<0.1
	medium	82	2.16	0.1
	high	98	1.76	0.2
Gastro, Gastro Standard Zoom 0, 30 x 30	low	90	0.44	<0.1
	medium	95	1.09	0.1
	high	99	1.75	0.2
Vascular, Vasc Standard Zoom 0, 30 x 30	low	71	1.28	<0.1
	medium	81	2.20	0.1
	high	97	1.83	0.2
Urology, Uro Standard Zoom 0, 30 x 30	low	71	1.40	<0.1
	medium	81	2.20	0.1
	high	97	1.83	0.2
Ortho/Trauma, Ortho Standard Zoom 0, 30 x 30	low	71	1.40	<0.1
	medium	81	2.20	0.1
	high	97	1.83	0.2
General, Standard Zoom 0, 30 x 30	low	71	1.41	<0.1
	medium	81	2.20	0.1
	high	97	1.83	0.2
General, Standard HL Zoom 0, 30 x 30	low	80	0.77	<0.1
	medium	82	2.17	0.1
	high	98	1.77	0.2
General, Standard HL (high dose) Zoom 0, 30 x 30	high dose	110	1.56	0.2
General, Pediatrics without scattered radiation grid Zoom 0, 30 x 30	low	85	0.17	<0.1
	medium	90	0.43	<0.1
	high	93	0.70	0.1

Important additional information

Dosimetric Information: valid from VA20H delivery state

Radioscopy, Detector size 20x20

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate at patient entrance reference point [mGy/min]
Cardiac, Card Standard Zoom 0, 20 x 20	low	80	1.5	33.3	cont	5,3
	medium	80	3.6	33.3	cont	12,7
	high	80	5.1	33.3	cont	17,5
Vascular, Vasc Standard Zoom 0, 20 x 20	low	62	24	10.4	10	4,3
	medium	66	24	14.5	15	10,4
	high	68	24	18.3	15	14,3
Urology, Uro Standard Zoom 0, 20 x 20	low	62	24	10.4	10	4,2
	medium	66	24	14.7	15	10,6
	high	68	24	18.3	15	14,3
Ortho/Trauma, Ortho Standard Zoom 0, 20 x 20	low	62	24	10.7	10	4,4
	medium	66	24	15	15	10,9
	high	68	24	18.3	15	14,4
Intervention, Vasc Interv Zoom 0, 20 x 20	low	64	24	14.3	7.5	4,7
	medium	66	24	15.4	15	11,2
	high	67	24	20.5	15	15,3
Intervention, Vasc Peripheral Interv Zoom 0, 20 x 20	low	64	24	13.7	7.5	4,5
	medium	66	24	15.4	15	11,2
	high	67	24	20.5	15	15,4
Intervention, Cardiac Interv Zoom 2, 20 x 20	low	80	14	7.1	15	4,9
	medium	80	3.6	33.3	cont	12,2
	high	80	5.1	33.3	cont	17,5
Intervention, Gastro Interv Zoom 0, 20 x 20	low	81	8.7	7.1	15	3,1
	medium	92	1.8	33.3	cont	8,9
	high	95	2.3	33.3	cont	12,2
General, Standard Zoom 0, 20 x 20	low	62	24	10.4	10	4,3
	medium	66	24	14.7	15	10,6
	high	68	24	18.6	15	14,6

Important additional information

General, Adipose Zoom 0, 20 x 20	low	63	24	10.7	10	4,4
	medium	66	24	14.7	15	10,7
	high	67	24	20.8	15	15,7
General, Pediatrics without scattered radiation grid Zoom 0, 20 x 20	low	69	6.3	7.1	10	0,9
	medium	76	7.1	7.1	15	2,1
	high	80	7.6	7.1	15	2,5
General, Extremity Zoom 0, 20 x 20	low	66	24	14.1	15	10,2
	medium	75	21	25.1	15	22,1
	high	84	18	25.1	15	26,0
General, Standard HL Zoom 0, 20 x 20	low	80	11	7.1	15	3,9
	medium	80	24	7.9	15	9,4
	high	80	24	16.3	15	19,8
Gastro, Standard Zoom 0, 20 x 20	low	86	1	33.3	cont	4,4
	medium	93	1.9	33.3	cont	9,2
	high	96	2.3	33.3	cont	12,3
General, Standard HL(high dose) Zoom 0, 20 x 20	high dose	110	8.9	33.3	cont	65.5

Radioscopy, Detector size 30x30

Program (with scattered radiation grid, Magnify 0	Dose level	Voltage [kV]	Current [mA]	Pulse width [ms]	Pulse rate [1/s]	Dose rate at patient entrance reference point [mGy/ min]
---	------------	-----------------	-----------------	------------------------	---------------------	--

Important additional information

Cardiac, Card Standard Zoom 0, 30 x 30	low	80	0.8	33.3	cont	2,8
	medium	80	2.1	33.3	cont	7,0
	high	80	3	33.3	cont	10,2
Vascular, Vasc Standard Zoom 0, 30 x 30	low	60	22	7.1	10	2,4
	medium	64	24	10.4	15	6,8
	high	65	24	13.1	15	9,1
Urology, Uro Standard Zoom 0, 30 x 30	low	60	22	7.1	10	2,4
	medium	64	24	10.4	15	6,8
	high	65	24	13.1	15	9,1
Ortho/Trauma, Ortho Standard Zoom 0, 30 x 30	low	60	22	7.1	10	2,4
	medium	64	24	10.4	15	6,8
	high	65	24	13.1	15	9,0
Intervention, Vasc Interv Zoom 0, 30 x 30	low	62	24	9.6	7.5	2,8
	medium	64	24	10.6	15	6,9
	high	65	24	14.1	15	9,4
Intervention, Vasc Peripheral Interv Zoom 0, 30 x 30	low	62	24	9.6	7.5	2,8
	medium	64	24	10.6	15	6,9
	high	65	24	14.1	15	9,4
Intervention, Cardiac Interv Zoom 2, 30 x 30	low	80	7.8	7.1	15	2,7
	medium	80	2.1	33.3	cont	7,0
	high	80	3	33.3	cont	10,2
Intervention, Gastro Interv Zoom 0, 30 x 30	low	75	6.9	7.1	15	1,9
	medium	89	1.3	33.3	cont	5,5
	high	91	1.7	33.3	cont	7,7
General, Standard Zoom 0, 30 x 30	low	60	22	7.1	10	2,4
	medium	64	24	10.6	15	6,9
	high	65	24	13.1	15	9,1
General, Adipose Zoom 0, 30 x 30	low	60	23	7.1	10	2,4
	medium	64	24	10.6	15	6,9
	high	65	24	14.1	15	9,4
General, Pediatrics without scattered radiation grid Zoom 0, 30 x 30	low	62	5.6	7.1	10	0,6
	medium	70	6.3	7.1	15	1,5
	high	73	6.8	7.1	15	1,8
General, Extremity Zoom 0, 30 x 30	low	64	24	8.1	15	5,3
	medium	66	24	21.1	15	15,0
	high	70	22	25.1	15	19,6

Important additional information

General, Standard HL Zoom 0, 30 x 30	low	78	7.3	7.1	15	2,3
	medium	80	15	7.1	15	5,2
	high	80	24	9.4	15	11,0
Gastro, Standard Zoom 0, 30 x 30	low	83	0.7	33.3	cont	2,5
	medium	89	1.3	33.3	cont	5,6
	high	91	1.7	33.3	cont	7,7
General, Standard HL(high dose) Zoom 0, 30 x 30	high dose	110	8.9	33.3	cont	67

Radiography, Detector size 20x20

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Charge [mAs]	Dose at patient entrance reference point [mGy]
Cardiac, Card Standard Zoom 0, 20 x 20	low	80	1.54	0.09
	medium	101	1.73	0.17
	high	110	1.59	0.19
Urology, Uro Standard Zoom 0, 20 x 20	low	73	2.52	0.11
	medium	100	1.78	0.18
	high	110	1.59	0.19
Ortho/Trauma, Ortho Standard Zoom 0, 20 x 20	low	73	2.25	0.11
	medium	100	1.78	0.18
	high	110	1.59	0.19
Intervention, Vasc Interv Zoom 0, 20 x 20	low	75	2.13	0.11
	medium	77	8.84	0.39
	high	77	11.60	0.59
Intervention, Vasc Peripheral Interv Zoom 0, 20 x 20	low	75	2.13	0.11
	medium	77	5.84	0.30
	high	77	11.60	0.59
Intervention, Cardiac Interv Zoom 0, 20 x 20	low	80	1.55	0.24
	medium	80	4.84	0.28
	high	80	9.41	0.34

Important additional information

Intervention, Gastro Interv Zoom 0, 20 x 20	low	93	0.77	0.17
	medium	96	2.15	0.20
	high	96	4.23	0.38
General, Standard Zoom 0, 20 x 20	low	73	2.52	0.11
	medium	100	1.78	0.18
	high	110	1.59	0.19
General, Adipose Zoom 0, 20 x 20	low	75	2.13	0.11
	medium	77	5.84	0.31
	high	77	11.60	0.59
General, Pediatrics without scattered radiation grid Zoom 0, 20 x 20	low	88	0.28	0.02
	medium	93	0.67	0.06
	high	96	1.17	0.11
General, Extremity Zoom 0, 20 x 20	low	87	1.97	0.14
	medium	110	1.59	0.19
	high	110	1.59	0.19
General, Standard HL Zoom 0, 20 x 20	low	80	1.59	0.09
	medium	101	1.73	0.18
	high	110	1.59	0.19
Gastro, Standard Zoom 0, 20 x 20	low	94	0.75	0.06
	medium	101	1.76	0.18
	high	110	1.59	0.19
Intervention, Vasc Interv (high dose) Zoom 0, 20 x 20	high dose	110	14	1.72

Radiography, Detector size 30x30

Program (with scattered radiation grid, Magnify 0)	Dose level	Voltage [kV]	Charge [mAs]	Dose at patient entrance reference point [mGy]
Cardiac, Card Standard Zoom 0, 30 x 30	low	80	1.03	0.03
	medium	87	2.03	0.15
	high	110	1.59	0.20

Important additional information


Urology, Uro Standard Zoom 0, 30 x 30	low	72	1.72	0.08
	medium	87	2.06	0.15
	high	110	1.59	0.20
Ortho/Trauma, Ortho Standard Zoom 0, 30 x 30	low	72	1.72	0.08
	medium	87	2.06	0.15
	high	110	1.59	0.20
Intervention, Vasc Interv Zoom 0, 30 x 30	low	73	1.55	0.07
	medium	75	4.30	0.21
	high	75	8.63	0.41
Intervention, Vasc Peripheral Interv Zoom 0, 30 x 30	low	73	1.55	0.07
	medium	75	4.30	0.21
	high	75	8.62	0.41
Intervention, Cardiac Interv Zoom 0, 30 x 30	low	80	1.07	0.06
	medium	80	3.23	0.19
	high	80	6.27	0.36
Intervention, Gastro Interv Zoom 0, 30 x 30	low	90	0.63	0.05
	medium	93	1.61	0.14
	high	93	3.14	0.26
General, Standard Zoom 0, 30 x 30	low	72	1.72	0.08
	medium	87	2.06	0.15
	high	110	1.59	0.20
General, Adipose Zoom 0, 30 x 30	low	73	1.55	0.07
	medium	75	4.30	0.21
	high	75	8.47	0.41
General, Pediatrics without scattered radia- tion grid Zoom 0, 30 x 30	low	86	0.21	0.01
	medium	91	0.49	0.04
	high	94	0.84	0.07
General, Extremity Zoom 0, 30 x 30	low	73	2.38	0.11
	medium	110	1.59	0.20
	high	110	1.59	0.20
General, Standard HL Zoom 0, 30 x 30	low	80	1.03	0.05
	medium	88	2.00	0.15
	high	110	1.59	0.20
Gastro, Standard Zoom 0, 30 x 30	low	91	0.54	0.05
	medium	97	1.30	0.12
	high	110	1.59	0.20

Important additional information

Intervention, Vasc Interv (high dose) Zoom 0, 30 x 30	high dose	110	14	1.75
---	-----------	-----	----	------

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
- ◆ Insert and secure the receiver module (1) in the DVI interface of your external monitor.
- ◆ Using a CAT cable (2), connect the receiver to the respective socket A (live monitor) or socket B (reference monitor) on the monitor trolley.
- ◆ Insert the 5 V connector of the power supply (3) into the receiver.
- ◆ Plug the corresponding power plug into the power outlet.
 - The connection to the external monitor is now operational.

HD Video manager	
Outputs	2 x CAT (RJ 45); DVI-D via active CAT/DVI adapter; 1280 x 1024 @ 60 Hz

Inputs	DVI-I (analog or digital) (2x), VGA with Adapter VGA to DVI-A, S-Video Mini-DIN and BNC, Composite, SD/HD-SDI/3g (2x), HDMI 2.0 (2x), Display Port 1.2; preprogrammed standard video timings up to (1920 x 1200, 50/60 Hz, interlaced/progressive) including PAL and NTSC DP and HDMI inputs video timing up to 4K@60 Hz
--------	--

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Risk of X-ray radiation without a stored image!

- ◆ In case of image failures stop exposures and switch off the system.

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Risk of electrical shock or damage to the system!

- ◆ Do not connect any additional devices to the UPS sockets.

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Radiation exposure without clinical benefit!

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CAUTION

Incorrect exposure parameters.

Risk of unnecessary radiation exposure due to incorrect exposure parameters!

- ◆ Ensure the exposure parameters are correct.
-

CAUTION

Storage rate is set incorrectly or being not aware of replay / image saving behavior

Images are lost or replay rate / image saving differs from the pulse rate during acquisition!

- ◆ Choose the storage rate with care.
 - ◆ Note that for 30 frames per second with 30x30 detector in MAG1 or with 20x20 detector in overview, image information is summarized in 15 frames for replay / storage only (no image information is lost).
-

CAUTION

The image orientation is reset when changing the application group or application.

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- ◆ Be aware that the image orientation is reset when changing the application group or application. All previous settings are lost.
-

CAUTION

Wrong calibration value is entered or length or angle information are not interpreted correctly.

Insufficient or incorrect basis for diagnostic or therapeutic decisions!

- ◆ Note that the measuring results of the functions Distance and Angle depend on various factors, e. g. the position of the calibration object. Note that the acquired images represent two-dimensional projections of real objects. Verify the results before diagnostic or therapeutic decisions are taken.

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CAUTION

Brake failure.

Crushing due to arbitrary movements!

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Spindle nut broken or worn.

Risk of crushing!

- ◆ Vertical lift is blocked if the mechanical safeguard was triggered. Contact Siemens Customer Service to replace the mechanical safeguard.

CAUTION

Audible warning for safety distance not functional.

Risk of crushing!

- ◆ Check for proper loudspeaker function before beginning the examination.
Note when booting the system whether an audible signal is heard.
-

CAUTION

Emergency stop not functional.

Risk of crushing and collision for persons and objects in the vicinity of moving parts!

- ◆ Check the EMERGENCY STOP function on a daily basis.
-

WARNING

Fluid or blood seeps into the system, use of inappropriate cleaning agents.

Risk of damage to the system, risk of infection or electrical hazard.

- ◆ Use sterile or sterilizable covers.
 - ◆ Do not let cleaning liquids seep into the openings of the system (e.g. air openings, gaps between covers).
 - ◆ Observe the cleaning and disinfection instructions.
 - ◆ Only use the recommended substances for cleaning and disinfection.
-

CAUTION

Under extreme operating conditions (continuous operation close to the maximum tube voltage), the X-ray tube can become very hot. This can shorten the service life of the system components.

Risk of burns and prolongation of the operation!

- ◆ Avoid touching the X-ray tube housing.
- ◆ Change your parameter settings, e.g. by reducing the dose level.
- ◆ We also recommend keeping a second system ready for continuous radiation lasting longer than 30 minutes.
- ◆ Please make sure that the system has cooled down before performing critical interventional applications.

WARNING

Fire in or near the system.

Injury to patient and personnel, and damage to device. Risk of gas poisoning due to burning plastic!

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WARNING

Due to the complexity of the system and related system component failures or disturbances in the mains power supply, e. g. line voltage is outside of the designated specification, it cannot be entirely excluded that system functions including X-ray imaging or other system functions get disturbed or unavailable.

Disturbance or unavailability of system functions. Interventions cannot be started or continued.

- ◆ Consider therefore the need to establish emergency procedures to manage situations related with the failures of the system, e. g. keeping a reserve system available.

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Risk of malfunction that could endanger the patient or product!

- ◆ Only use Siemens-authorized software and hardware components.
 - ◆ Only perform repairs, or have them performed by an authorized third party, after having received our express written approval.
-

WARNING

Cable damage.

Risk of electrical shock!

- ◆ Check the power cable. Do not use the device if the power cable is damaged.
-

CAUTION

Image intensifier / flat detector cover broken

Skin or eye irritation, irritation of the breathing passages, toxic or allergic reactions caused by contact with hazardous substances!

- ◆ If the image intensifier / flat detector cover is broken, do not touch the detector, do not use the system at all.
 - ◆ Please call Siemens service.
-

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Risk of crushing!

- ◆ Only products/components approved by Siemens may be installed.
-

WARNING

Interference.

Impact on the patient's electronic life-sustaining systems!

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WARNING

The system is not designed for operation in explosion-endangered areas. It does not fulfill the requirements for AP/APG classification.

Explosion hazard!

- ◆ The system may may not be used in locations where explosive atmospheres may occur.

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Misinterpretation of image information or of indications of the system, leading to inappropriate decisions in the clinical workflow, incorrect diagnosis or treatment, delay of the clinical procedure, or exposure to radiation without clinical benefit. Only persons with the required knowledge and expertise who have undergone the appropriate user training may work with the system.

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CAUTION

Unintentional activation of movement

Risk of injury and crushing!

- ◆ Angulating the C-arm > 90° may cause the cable to collide with the operating unit on the horizontal support arm and accidentally operate movements.
-

CAUTION

Processing images during a lithotripsy procedure.

Risk of incorrect decision during intervention, unnecessary radiation exposure and unnecessary stress on tissue and/or organs by shock-wave application!

- ◆ Do not process images in lithotripsy mode.
-

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- ◆ Make sure no instruments, conductive material, or fluids penetrate into the openings/joints between the control unit and its housing.
-

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Damage to the flat detector!

- ◆ Make sure that the detector is always covered by a disposable cover during the examination.
-

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Risk of injury!

- ◆ Use only approved accessories and combined devices.

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Risk of injury!

- ◆ Note that stability in accordance with DIN EN/UL standards is guaranteed only if the external system (navigation) on the detector does not exceed the calculated maximum weight.

CAUTION

Electromagnetic interference of the system.

Degradation of the image quality or accuracy of the radiation parameters!

- ◆ Remove the source of interference (for instance, a mobile phone) from the area.
- ◆ Shut down the system and turn it on again.
- ◆ Pay attention to error messages during restart.

Important additional information

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Addendum to the User Manual for Russia

Mobile X-ray C-arm Unit for radioscopy “Cios Fusion” with accessories

1. The name of the Medical Device.

Mobile X-ray C-arm Unit for radiscopy "Cios Fusion" with accessories

2. Manufacturing information of the Medical Device.

Legal manufacturer and Developer:

Siemens Healthcare GmbH
Henkestr. 127,
91052 Erlangen, Germany

Manufacturing sites:

Siemens Healthcare GmbH
Advanced Therapies
Röntgenstr. 19-21, 95478
Kemnath, Germany

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubininskaja street 96,
115093 Moscow, Russia
Tel.: +7 495 737 1252
www.healthcare.siemens.ru

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3.3. User profile

The system can be used properly only by suitably qualified medical personnel. This requires operating personnel to be familiar with the Manual. This manual must be studied in detail prior to starting up the system.

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As the manufacturer of the system, SIEMENS is not responsible for any impact on its safety, reliability and performance if it is not used in accordance with the Manual.

3.4. Residual risk

The Cios Fusion is a complex medical device, its functionality requires the seamless and synchronized operation of a high number of hardware and software components, as well as appropriate operating conditions.

A safe and effective operation of the Cios Fusion requires that it is operated by persons with the necessary specialist knowledge and appropriate knowledge about the system, its intended use, its functions, the conditions and limitations for its usage, periodic maintenance and routine checks which need to be performed, the potential error situations, as well as the corresponding recovery possibilities.

System functions may be disturbed or completely lost at any time due to component errors, damages like ingress of liquids or collisions, or due to inappropriate operating conditions e.g. with regards to mains power supply, excessive use of the system at maximum load for a prolonged period of time or inappropriate maintenance or service. User or operator behavior like using the system without appropriate care, using the system in combination with devices, accessories, other equipment or pieces of software in a way which is not approved by the manufacturer or cleaning or disinfecting the system not according to the provisions given in the instructions for use may also lead to damages or functional disturbances. Certain materials in the X-ray beam like patient tables may impair the image quality and may lead to a higher radiation exposure.

The Cios Fusion is designed to comply with regulatory requirements regarding electromagnetic compatibility. It is protected against the electromagnetic influence of external devices and does not disturb the operation of these. However, individual external devices with emissions beyond the allowed levels or even the combined emission of multiple, individually compliant devices may lead to the influence levels considered in designing and testing the Cios Fusion being exceeded, and eventually, to functional disturbances. External devices, which are not sufficiently stable against electromagnetic influence may be disturbed by the emissions of the Cios Fusion as well.

When releasing system movements, the operator's care is required in order to prevent crushing, collisions with other equipment or collisions with and thus movement of body-penetrating objects like catheters or needles.

The Cios Fusion manages radiation exposure in order to minimize the radiation dose for the patient, user and staff, while providing the necessary level of imaging performance. The applied dose is measured, tracked and displayed. Warnings are given under defined conditions.

In combination with the appropriate care in the clinical workflow and informed decisions considering also the clinical benefit of radiation, the above mechanisms should prevent deterministic dose effects to the patient (erythema, cataracts, permanent epilation and delayed skin necrosis). Nevertheless, such

injuries cannot be completely excluded. In rare cases, the dose management mechanisms may be not be operating correctly. These issues are detectable for the user especially within a Quality Assurance framework. Some clinical applications using the Cios Fusion involve injecting contrast medium into the patient's cardiovascular system or other body orifices. Functional disturbances in the Cios Fusion may affect the imaging workflow resulting in the additional administration of contrast medium. This may lead to risks and side-effects in accordance with the information provided by the manufacturer.

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3.6. Undesireable side effects:

Deterministic and stochastic effects of exposure to ionizing radiation (X-rays) are part of the physician's medical education. The necessity of an X-ray examination in relation to these risks must be considered by the physician when prescribing or performing such a type of examination.

For radiation emitting products, known deterministic effects of exposure to X-rays are erythema, cataracts, permanent epilation and delayed skin necrosis.

For radiation emitting products, known stochastic effects of exposure to X-rays are increased risks of developing cancer and hereditary diseases.

These stochastic dose effects have an increased probability of occurrence with increased dose. The related risk can be minimized by means of keeping the X-ray exposure to the necessary minimum and by means of careful and systematic application of protection measures for the operator and the clinical staff.

Examinations performed on patients especially sensitive to ionizing radiation (e.g. infants, pregnant women, individuals with certain genetic predispositions and/or diseases characterized by enhanced radiation sensitivity) or to contrast medium (e.g. chronic kidney diseases) need to be planned and performed with special care to minimize the exposure.

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Work environment:

Temperature range: +15 °C to +35 °C

Rel. humidity: 15% to 75% (without condensation)

Barometric pressure: 700 hPa to 1060 hPa

4. Classification of the Medical Device:

4.1. Type of the medical device in accordance with the nomenclature for classification:

Classification according to Annex IX, MDD

Class	I Ib
Classification in accordance with Russia Nomenclature classification:	

209270

4.2. Class of potential Risk of the medical device application:

The Medical Device is intended for transmission of energy in the form of ionizing radiation, it belong to the **Class 2b**.

4.3. Class of the Medical Device in accordance with IEC 60601-1:

Protection against electric shock: Class 1, no applied part according to IEC 60601-1

Operation mode: Pulsed and continuous radioscopy

4.4. Type on contact with the human body:

The Medical Device is non-sterile, nonexpendable.

Device surfaces that can come into contact with patients are evaluated and proven to be biocompatible according to ISO 10993:2009.

4.5. IPX classification

X-ray tube assembly (Single-tank)	IPX0 - No protection
Flat detector	IPX0 - No protection
C-arm	IPX0 - No protection
Remote control unit (optional)	IP65 dust tight, protected against water jets
Standard footswitch	IPX8 protected against immersion beyond 1 m
Multifunctional footswitch	IPX8 protected against immersion beyond 1 m
Hand switch	IP54 dust protected, protected against splashing of water

4.6. Laser positioner:

Laser light Class: **1** (according to IEC 60825:1:2007).

4.7. Software classification:

The Medical Device's Software safety Class according to IEC 62304:2006 is **Class B**.

4.7.1. Software version

Software version	Software release version: VA20x. Date of release VA20x: July 2015 Date of release VA20H: June 2021 Full software version: VA20C/VA20D/VA20E/VA20G/VA20H and later.
------------------	---

Additionally, for Cios Fusion systems with the VA20x software version previously registered in the Russian Federation (according to the registration certificate RZN 2018/7914 dated 10/12/2018), it is possible to upgrade to the following software patches:

VA20A	VA20E
VA20B	VA20F
VA20C	VA20G
VA20D	VA20H

All information related to the VA2x software version is indicated in the registration dossier RD-22112/15289 dated 11/05/2018 (registration certificate to the RZN 2018/7914 certificate dated 10/12/2018). No additional technical parameters and characteristics are required for the above software patches (see information on software version naming below).

The versioning of software is handled according to following definitions and abbreviations.

The Software Version serves to distinguish between different software statuses:

V means released version

A means uppercase letters A to Z without O, I and without language-specific characters (e.g.: Ä, Ď, E, Ö),

N means only numbers 0 (zero) to 9 (nine)

VANNA criteria:

Main Version

The main version is changed to the next letter, if software is not backwards compatible or if a sellable

option (possibly configurable) was changed. The change of the main version has to lead into a change of material number. Otherwise the material number of the software is not changed.

Sub Version

The sub version is changed to the next number if new functionality is implemented.

Patch Level

The patch level is changed to the next letter due to correction of errors.


4.8. *SanPiN 2.1.3684-21 class*



Class A

5. Technical Data: Specifications and Characteristics.

Accuracy: $\pm 5\%$ (unless otherwise indicated)

5.1. Dimensions of the system (cm)

5.2. Technical Data

General data

Power requirements	100 V, 110 V, 120 V, 127 V, 200 V, 220 V, 230 V, 240 V $\pm$ 10%; 50/60 Hz ($\pm$ 1 Hz) An on-site 16 A/characteristic C overcurrent release (trip breaker) is recommended by DIN VDE 0100-710 for supply network voltages of 230 V/240 V. Country-specific requirements apply for the fuse values to be used for other supply network voltages.
Unit fuse protection (internal)	100 V to 127 V 20 A slow-blow fuse 200 V to 240 V 15 A slow-blow fuse
Nominal rating	20 A to 127 V~, 15 A from 200 V corresponds to the nominal value of the slow-blow fuse in the power input of the product
Internal line impedance	Ri max. 0.3 Ohm for 100 - 127 V Ri max. 0.8 Ohm for 200 - 240 V
Power consumption	Maximum 2.6 kW Standby 475 W (for 230 V)
Continuous	14.5 A (100 V) / 7 A (230 V)
Short term (max. 30 s 1000 W Fluoro)	25 A (100 V) / 12 A (230 V)
Short term (max. ~ 50 ms DR single shot)	29 A (100 V) / 14 A (230 V)

Weight of the main components

Monitor trolley (including 2 monitors, UPS)	230 kg
---	--------

C-arm system without accessories	272 kg
----------------------------------	--------

Current/voltage values

Voltage (V AC)	Long-term current consumption 50 Hz/60 Hz (A)	Short-term current consumption 50 Hz/60 Hz (A)
100	18	29
110	17	27
120	16	25
127	16	24
200	10	16
220	9	14
230	9	14
240	9	14

X-ray generator

High-frequency generator	2.3 kW (non-IEC)
Rated power (IEC 60613)	1.4 kW (100 kV, 14 mA, 100 ms)
Single image pulse power (max)	1.9 kW (80 kV, 23.7 mA, 33.5 ms)
Average power in series operation (max)	1000 W (71 kV, 14.1 mA, 33.3 ms, 30 P/s)
Inverter control frequency	15 kHz to 30 kHz
kV range	40 kV to 110 kV
Single image	0,2 mA to 25 mA 0.021 mAs to 4 mAs 7 ms to 160 ms
Continuous fluoroscopy	0,2 mA to 15 mA
Fluoroscopy	3 mA to 25 mA 0.5 P/s to 15 P/s (30 P/s for continuous fluoroscopy) 7 ms to 40 ms pulse width
Single tank heat storage capacity	888,000 J/1,200,000 HU

Data according to IEC 60601-2-54:

Radiography	
Rated tube voltage	110 kV
Highest X-ray tube current	25 mA
Highest power (nominal power)	100 kV, 14 mA, 100 ms
Smallest value for current-time product	0.02 mAs
Data according to IEC 60601-1-3 Accuracy of loading parameters	Tube voltage (kV): $\pm 5\%$ Tube current (mA): $\pm 15\%$ (minimum ± 0.1 mA) Current-time product (mAs): $\pm 10\%$ (minimum ± 0.2 mAs) Pulse duration (ms): $\pm 10\%$ (minimum ± 1 ms)

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Model	Single tank with single-focus stationary anode tube
Tube material	Tungsten
Nominal focus value	0.6
Tube rated voltage (IEC 613)	110 kV
Single tank rated voltage	110 kV
Anode heat storage capacity (IEC 613)	45,200 J (61,100 HU)
Anode heat dissipation	27,600 J/min (37,300 HU/min)
Anode angle	9°
Inherent filtration tube (IEC 601)	2.8 mm Al with 75 kVp/0.1 mm Cu
Total filtration system (IEC 601)	3.0 mm Al with 75 kVp/0.1 mm Cu
Continuous heat dissipation	80 W
Continuous input power	170 W
Maximum irradiation time	20 min at 600 W / 60 min at 170 W

Flat detector

Cios Fusion utilizes Varian PaxScan detectors. The solid state detector converts the X-rays to visible photons by a scintillator. The visible photons are absorbed by pin photodiodes and converted to electron-hole pairs, which collect on the capacitance of the photodiodes. The pin photodiodes are discharged when the pixel's TFT is turned ON. The charge is collected by an integrating amplifier and converted to a voltage. The signal voltage has a programmable gain applied, depending upon operating conditions. The output voltage is converted to digital data by an analog-to digital converter.

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	dose (fluoro)							
	Saturation dose (single image)	μGy	29			29		
	Noise equivalent dose	nGy/ binning		3.5/2x2			3.5/2x2 ¹⁰⁰	
	Dynamic range 12)	dB		94	96		94	96
Image quality performance 13)								
	Sensitivity (low-dose fluoro)	LSB/nGy	30	36		17	20	
	Electronic Noise (lowdose fluoro)	nGy/ binning			0.75/2x2			1.5/1x1
	Residual Lag	% in 1 st image			5			5
DQE @ 40nGy, RQA5 (30x30 cm2, 2x2)								
	at 0.05 lp/mm	%		76			n/a	
	at 1.0 lp/mm	%		43			n/a	
	at 1.29 lp/mm (Nyquist)	%		25			n/a	
DQE @ 80 nGy, RQA5 (20x20 cm2, 1x1)								
	at 0.05 lp/mm	%		76			76	
	at 1.0 lp/mm	%		55			55	
	at 2.0 lp/mm (Nyquist)	%		30			30	
	at 2.58 lp/mm (Nyquist)	%		14			14	
MTF (1x1 binning)								
	at 1.0 lp/mm	%		55			55	
	at 2.0 lp/mm (Nyquist)	%		25			25	
	at 2.58 lp/mm (Nyquist)	%		16			16	
Main functionalities								
	Max frame rate for full field of view ¹⁴	fps/binning	15/1x1		30/2x2			30/1x1
	Max frame rate for limited field of view ¹⁵	fps/binning			30/1x1			n/a

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Scattered radiation grid on flat detector

Varian PaxScan 2020XV	Pb 17 / 70 fo 100
Varian PaxScan 3030XV	Pb 17 / 70 fo 100

C-arm

Orbital movement	150° (-40° to +110°)
Angulation	± 225°
Horizontal movement	20 cm (7.9")
Immersion depth	73 cm (28.7")
Swivel range	± 12°
Vertical movement	43 cm (16.9"), motorized
SID	102 cm (40.2")
Single tank image receptor distance	82 cm (32.3")

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Monitors

TFT high-brightness color display

Screen diagonal	19" (48 cm)
Image display	1280 x 1024 (pixels)
Brightness (DICOM calibrated)	400 cd/m2
Maximum brightness, typical	650 cd/m2
Horizontal/vertical viewing angle	178°/178°
Contrast ratio, typical	900 : 1
Backlight LED technology	

TFT Premium high-brightness color display

Anti-reflexion coated cover	Glass panel for higher-contrast image quality and easier cleaning
Brightness Uniformity Auto Adjustment	Brightness uniformity is automatically adjusted for a more balanced and uniform image impression
Automatic Brightness Adjustment Brightness is automatically adjusted based on ambient lightAutomatic DICOM Correction	DICOM curves are automatically corrected based on ambient light (This feature greatly speeds up verification of conformity to DIN 6868-157 in affected countries)
Full b/w emulation mode	
Protective collision cover	
Monitor Interface	Input: 1 x DisplayPort
Screen diagonal	19" (48 cm)
Image display	1280 x 1024 (pixels)
Brightness (DICOM calibrated)	500 cd/m2
Maximum brightness, typical	1000 cd/m2
Horizontal/vertical viewing angle	170° / 170°
Contrast ratio, typical	1000 : 1
Response time, typical	25 ms
Backlight LED technology	

Monitor columns

Flex display column (Standard)	Display positioning with rotation angle -30° to 210°; defined lock-in positions at 0° and 180°
--------------------------------	--

Flex Plus display column (optional)	Display positioning with rotation angle -30° to 210°; defined lock-in positions at 0°, 90° and 180°, motorized height adjustment; fold function to protect the displays in park position and during transport
-------------------------------------	---

Dose measuring chamber (DAP meter)

Technology	Semiconductor chamber
Active area	Rectangle - 60 mm x 60 mm
DAP resolution	0.01 μGym^2 equivalent to 0.01 cGycm ²
Maximum measurable DAP	3000.00 μGym^2 equivalent to 3000.00 cGycm ²
Reproducibility	< 1% under constant pressure and temperature
Energy range	40 - 150 kV $\pm$ 6%
Attenuation equivalent (Inherent filtration)	< 0.30 mm Al equivalent
Dose rate linearity	Better than $\pm$ 2% for dose area product rate $\geq$ 1 $\mu\text{Gym}^2/\text{s}$, otherwise better than $\pm$ 5%
Measuring readiness	20 seconds after power on
Work environment	Temperature range: +10 °C to +70 °C Rel. humidity: 15% to 75% (without condensation) Barometric pressure: 700 hPa to 1060 hPa
Precision of the dose measuring chamber (DAP meter)	< 25% (deviation between displayed and measured dose area product)

Automatic dose rate control (ADR)

Precision of automatic dose rate control (ADR)	$\pm$ 25% (value depends on measurement setup)
--	--

Detector laser light localizer

Laser class	Class 1 (IEC 60825-1:2007)
Laser type	Semiconductor laser (laser diode)
Wave length	650 nm
Color	Red
Effective power	$\leq$ 3 mW

Single-tank laser light localizer

Laser class	Class 1 (IEC 60825-1:2007)
Laser type	Semiconductor laser (laser diode)
Wave length	650 nm
Color	Red
Effective power	≤ 3 mW

HD video manager

Outputs	2 x CAT (RJ 45); DVI-D via active CAT/DVI adapter; 1280 x 1024 @ 60 Hz
Inputs	DVI-I (analog or digital) (2x), VGA with Adapter VGA to DVI-A, S-Video Mini-DIN and BNC, Composite, SD/HD-SDI/3g-SDI (2x), HDMI 2.0 (2x), DisplayPort 1.2; preprogrammed standard video timings up to (1920 x 1200, 50/60Hz, interlaced/progressive) including PAL and NTSC DP and HDMI inputs provide video timing up to 4K@60Hz

DVI video splitter


Outputs	2 single-link DVI-D; female connectors; up to 1920 x 1200 @ 60 Hz or 1080 p
Conditions for use	The displays or monitors used must be compatible with the SXGA standard 1280 x 1024 at 60 Hz
	No galvanic separation, therefore the connected monitors must comply with IEC 60601-1-1

SCALANCE W 721 WLAN access point

Data transfer	
Ethernet transfer rate	BASE-T max. 100 MBit/s
Wireless transfer rate	max. 150 Mb/s
Radiated power	typically: 100 mW, max. 200 mW

Supported wireless standards	IEEE 802.11a IEEE 802.11b IEEE 802.11e IEEE 802.11g IEEE 802.11h IEEE 802.11i IEEE 802.11n
Supported standards for power supply	IEEE 802.3at Type 1 (802.3af) and 802.3at Type 2 (power supply via Ethernet)

6. Scheme of the Production line.

Quality check

- System check by Designated Q-Person SCM (includes a check of functionality, a check of completeness of customer order, test protocols, customer documentation and system labeling)

System disassembly

- Move the system to transport position for applicable components
- Disconnect system cabling
- Assemble transport racks
- Prepare and provide system components for packing

Hand over to logistic provider


7. The list of components and accessories.

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



камера для измерения дозы облучения	Dose measure. chamber	System-integrated dose measuring chamber for displaying the dose area product or air kerma value. The cumulative dose area product is displayed for the current patient and saved under the patient data. The cumulated dose is automatically transferred to a radiation summary report and can be retrieved at any time. For each patient a cumulative value is saved in the patient database.	
документация пользователя на русском языке	User documentation	The manual is intended for health professionals with training in the field of radiology and aims to familiarize them with the system components and functions. It is not intended to train them in radiology or make any clinical diagnoses.	
принадлежность	Accessories:		
Лазерный оптический центратор детектора	Int. detector laser light localizer	Integrated laser aimer on the detector which transmits a crosshair laser. For precise, radiation-free positioning with tube in undertable position	



2. Моноблочный лазерный оптический центратор	Int. single tank laser light local	Integrated laser light localizer in single tank which emits a crosshair laser. For precise, radiation-free positioning in overtable position.	
3. Программный режим Памяти положения	Position memory	Position memory for storing and retrieving the orbital/ and angular position easily. The position memory provides an additional key on the touch control panel which enables the orbital and /angular position to be stored and retrieved easily. The clear display of the angular and orbital angles with orientation arrows on the control panel allows the stored value to be retrieved easily.	
4. Разветвитель видеосигнала DVI	DVI video splitter 2x	Connection for an external live monitor (monitor A) and an external reference monitor (monitor B) via DVI connection.	


Диспетчер видеонаблюдения HD (SmartView)	SmartView - HD VideoManager	HD video manager for emitting video signals, e.g. on external monitors, or for feeding in video signals for display on the reference monitor (monitor B), e.g. endoscopy. Interface for display of external video signals such as endoscopy on the right-hand TFT color display on the unit. Supports HD (High Definition) endoscopy.	
Встраиваемая акустическая стереосистема	EMotion sound system	EMotion integrated sound system with interface for operating external (MP3) audio equipment via Aux-in jack (Aux input) including amplifier with digital sound processor (DSP) and 2-way stereo loudspeaker system for high-fidelity sound.	
Пульт дистанционного управления с держателем	Remote control unit incl. holder	Touch-based remote control unit for operating the C-arm from within the sterile work area. C-arm functions that previously could only be controlled from within the non-sterile area, are now accessible from within the sterile area via the remote control unit.	

8. Дополнительный держатель для пульта дистанционного управления	Addit. remote contr. unit holder	An additional table holder requires an existing, complete remote control unit.	
9. Программный режим 2D измерения (для измерения углов и расстояний)	2D measurement function	Measurement of angles and distances. 2D measuring function with integrated calibration for determining the scale for measuring lengths and distances in the X-ray image. When measuring angles, the sides of an angle can be changed independently of one another. Complementary angles (difference of 180° - measured angle) can be displayed by changing the angle direction.	

Програм- мный режим для обсле́дова- ния сосудов	Vascular Software	Acquisition rate can be selected up to 30 f/s, storage rate can be configured. Automatic request for contrast medium injection following generation of the mask image. Display of dynamic changes during the contrast medium flow. Further processing functions: - Pixelshift - Remask	
Програм- мный режим для обсле́дова- ния сосудов Premium	Vascular Software Premium	Acquisition rate can be selected up to 30 f/s, storage rate can be configured. Automatic request for contrast medium injection following generation of the mask image. Display of dynamic changes during the contrast medium flow. Further processing functions: - Pixelshift - Remask - Draw AAA Line	


12. Программный режим активного наложения графики (построение линий и отображение на прямом изображении)	Live Overlay	Graphical	Optional software for drawing lines in images during the examination. The lines drawn in an ongoing examination can be transferred to the subsequent live images.	
13. Программный режим Количественное представление стеноза	Stenosis quantification		Quantification program for determination of degree of stenosis Calculation of the degree of vessel stenoses. Automatic vessel analysis and display of the minimum value in relation to a defined reference value. Geometrical and densitometrical determination of the value.	
14. Разъем для инжектора контрастного вещества	Injector trigger		Unidirectional trigger output for connecting a contrast agent injector.	



Програм- мный режим для записи рентгеноскоп ических последовател ьностей на CD/DVD- диске в формате MPEG2	CD/DVD Recorder_Subs	Software license for direct output of fluoroscopy and image series on DVD recorder in MPEG2 format.	
Програм- мный режим DICOM Send (Отправка DICOM)	DICOM Send/Storage Commitment (StC)	Digital, unidirectional image transfer of single images or complete folders to a network in DICOM format. Feedback from the image archive (Storage Commitment).	



17. Програм- мный режим DICOM Print (Печать DICOM)	DICOM Print	For sending and printing of images by means of a virtual filmsheet to a DICOM laser camera or printer. Provision of DICOM Print service for connection to a laser camera or a network printer (postscript- capable).	
18. Програм- мный режим DICOM Query/ Retrieve (Запрос /получение DICOM)	DICOM Query / Retrieve	Retrieval of archived images (DICOM Query/Retrieve - Service as SCU) from a digital archive or a workstation. The images must have a DICOM XRF (X-Ray Radiofluorography)/XA (X-Ray Angiography), CR (Computed Radiography), or SC (Secondary Capture) format.	

Програм- мный режим DICOM Worklist (Рабочий список DICOM)	DICOM Worklist / MPPS	Import of patient/examination data from an external RIS/HIS patient management system with DICOM MWL (Modality Worklist) as well as feedback on the examination status with DICOM MPPS (Modality Performed Procedure Step).	
Програм- мный режим DICOM Advanced	DICOM Advanced	Advanced DICOM is a package comprised of the following single options: - DICOM Storage Commitment - DICOM Print - DICOM Modality Worklist - DICOM Query/ Retrieve	
Цифровой интерфейс для передачи без потерь данных двухмерного ображе- ния в навигацион- ную систему - NaviLink 2D	NaviLink 2D	Integrated navigation interface for connecting navigation systems with digital DICOM 3.0 image transfer to the navigation system 2D Image transfer	

22. Принтер для вывода изображений на бумагу/плёнку	Sony Printer UP-971AD Sony Printer UP-991AD, 971 AD Sony Printer UP-D898MD	Thermal printer for gray-scale printing	
23. Устройство для установки принтера	Printer installation kit	Installation kit for connection of the printer to the Cios Fusion	
24. Программный режим электронной блокировки излучения	Radiation locking code	Password-protected radiation block. Radiation block with configurable password / PIN to avoid accidental radiation release.	

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
Main systems label:

Artwork of the MD labeling in Russian

Additionally, for all executions (variants) of medical device that are delivered on the territory of Russian Federation, a Russian label is applied with the following content:

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
8.1. Location of labels.


C-arm system

The labels shown below are attached permanently to the following subassemblies.

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












Multifunctional footswitch


(8) Multifunctional footswitch identification label

8.3. C-arm system control unit and remote control unit

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


8.4. Flat detector, grid, frame, laser light localizer

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

(21) xxxxx SERIAL

IVK

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured:
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

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(21) xxxxx SERIAL

IVK

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured:
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

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Made in Germany

Identification label


Warning label (hot surface)
The label is attached to the front, left, and right side of the single tank.

X-ray tube


Identification label

Approval-related label, filter label
Collimator

Identification label

Approval-related label


Generator

Identification label

8.6. Single tank laser and dose measuring chamber

(1) Single tank laser


This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture.

Manufactured:

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

(2) Dose measuring chamber

This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture.

Manufactured:

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

8.7. Monitor trolley

component name

(01) xxxxxxxxxxxxxxxxxxxx	GTIN
(240) xxxxxxxx	MODEL
(422) xxx	ORIGIN YY
(21) xxxxx	SERIAL

IVK

This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture.

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S.N: AN001152000000

100-240V ~ 50/60Hz 0.7 0.3A

DISCONNECT SUPPLY BEFORE SERVICING

One device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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WARNING!

THIS X-RAY UNIT MAY BE DANGEROUS TO PATIENT AND OPERATOR UNLESS SAFE EXPOSURE FACTORS, AND OPERATING INSTRUCTIONS AND MAINTENANCE SCHEDULES ARE OBSERVED.


(4) Warning label

			Long Time Current	Moment Current
Volt	50 Hz	60 Hz	A	A
100 ~	☐	☐	18	29
110 ~	☐	☐	17	27
120 ~	☐	☐	16	25
127 ~	☐	☐	16	24
200 ~	☐	☐	10	16
220 ~	☐	☐	9	14
230 ~	☐	☐	9	14
240 ~	☐	☐	9	14

(1) Supplemental label country-specific connection data

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(4) Approval label (not for China)

(5) Warning label read Operator Manual, watch inclines

9. The list of International Standards and Directives applied to the Medical Device.

Standard organization and reference number and date	Title of Standard
IEC 60601-1 3rd	Medical Electrical Equipment - Part 1: General Requirements for Safety, 2005
IEC 60601-1-1:2000	Medical electrical equipment – Part 1-1: General requirements for safety - Collateral standard: Safety requirements for medical electrical systems
IEC 60601-1-2 3rd 2007	Medical Electrical Equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
IEC 60601-1-3: 3rd 2008	Medical electrical equipment – Part 1-3: General requirements for basic safety and essential performance - Collateral Standard: Radiation protection in diagnostic X-ray equipment
IEC 62366: 2007	Medical devices - Application of usability engineering to medical devices
ISO 14971: 2012	Medical devices - Application of risk management to medical devices
IEC 62304 Ed. 1.0 2006	Medical device software - Software life cycle processes
IEC 60601-2-7 (1998)	Medical electrical equipment – Part 2-7: Particular requirements for the safety of high-voltage generators of diagnostic X-ray generators
IEC 60601-2-28 Edition 2.0 2010-03	Medical electrical equipment – Part 2-28: Particular requirements for the basic safety and essential performance of X-ray tube assemblies for medical diagnosis
IEC 60601-2-32: 1994	Medical electrical equipment - Part 2: Particular requirements for the safety of associated equipment of X-ray equipment
IEC 60601-2-54 2009	Medical electrical equipment – Part 2-54: Particular requirements for the basic safety and essential performance of X-ray equipment for radiography and radioscopy
IEC/PAS 61910-1 2007	IEC PAS 61910-1 First Edition 2007-07, Medical electrical equipment - Radiation dose documentation – Part 1: Equipment for Radiography and Radioscopy
NEMA PS 3.1 - 3.20 (2011)	Digital Imaging and Communications in Medicine (DICOM) Set
IEC 60825-1:2007	Safety of laser products - Part 1: Equipment classification, and requirements

IEC 60601-2-43 - Ed. 2.0 2010-03	Medical electrical equipment – Part 2-43: Particular requirements for the safety and essential performance of X-ray equipment for interventional procedures
-------------------------------------	---

10. Sterility.

The Medical device is non-sterile (except for sterile covers).

11. Disposal.

	The use of this symbol indicates that this product should not be treated as household waste. By ensuring that this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste-handling of this product. For more detailed information concerning the return and recycling of this product, please consult the supplier from whom you purchased the product
	Concerning Batteries (Directive 2006/66/EC) The directive 2006/66/EC requires separate collection and appropriate disposal of spent batteries. This product contains batteries that are not intended to be replaced by the user. Replacement of those batteries will usually be done during regular maintenance or service by service staff who can also arrange proper disposal.

If you want to remove the product from service, take into consideration that public legal directives may contain special regulations regarding disposal of this equipment. In order to ensure that these legal regulations are complied with and to avoid potential environmental hazards which may be caused by the disposal of your system, please consult Siemens Customer Service.

- Batteries and packaging material must be disposed of in an environmentally safe manner according to national regulations.
- The sterile single-use covers must be disposed of in accordance with national regulations or the rules of the hospital.
- For further information about the disposal of the product, please refer to our service documents.

12. Packaging.

The packaging consists of wood, nails, cardboard, paper and plastic. The packing list indicating names, quantities, packaging date, and packer reference number is inserted in the shipping box.

The following items appear on the label affixed to each package:

- Manufacturer trademark;
- Device name;
- Manufacturing date.

The following items appear on the label made according to the manufacturer's working drawings and affixed to each multiple shipping box:

- Medical device manufacturer name;
- Net and gross weight;
- Box dimensions

Transport marking. The shipping box is marked with handling symbols which have the following meanings: "This side up", "Keep away from moisture", "Fragile, handle with care!".

Packaging labeling

	"Fragile" Indicates that medical device can get broken or damaged if not handled carefully.
	"This side up"
	"Keep away from moisture" Indicates that the medical device has to be protected from moisture.

C-arm packaging process

- Packaging and transportation of C-arm is performed after it's positioning on transport platform.
- Than special support unit is mounted which is secured by brackets. Also other parts of C-arm are secured by brackets.
- Bracing cables and mounting brackets secure box with printer and documentation.
- The packing list is inserted in the shipping box. Package is covered by protective layer and secured by yellow ramps.

Monitor trolley packaging process.

- After all wheel lock ups unlock monitor trolley is rolled on transportation platform and secured by special foot blocks.
- Wheels are turned along the line and secure them with brakes.

- Foot blocks positioned at right and left sides provides secure of medical device on transportation platform.
- Secured monitor and keyboard are covered with protective layer which protects them from damage during transportation.
- The packaging process is completed with securing side panel boards and labeled cover up (see beginning p.9).

Real photos of the packaging:

Dimension of packaging of Cios Fusion (Length x width x height) : 2470 x 1150 x 1990 mm.
Cios Fusion with packaging total weight: 648kg.

13.Storage requirements and conditions.

Storage temperature range: -20 °C - +70 °C
Rel. humidity: 10% to 95%, non-condensing
Barometric pressure: 700 hPa to 1060 hPa

14.Service life.

Our medical engineering products are designed for operation under average conditions in accordance with the Operator Manual for a useful life of **10** (ten) years. If the products are operated for longer than this time, additional checks and possibly repairs beyond the usual maintenance procedures may be necessary in order to ensure the functional integrity and operating safety of these products.
Please contact Siemens about these measures early enough!

15. Warranty.

Siemens guarantees that the Equipment has no defects caused by installation and, as of the date of risk transfer, is compliant with the approved specifications and has no defects in materials and production within 14 (fourteen) months as of the Equipment delivery date. Warranty terms and conditions for components made by the third parties and supplied along with the Equipment shall be detailed in the respective Order.

Warranty period for operable Equipment components supplied to replace the defected ones shall be 6 (six) months as of their delivery date. Total warranty period for defective components repaired or installed by Siemens shall not exceed 14 months as of the Equipment delivery date.

Defective Equipment (its part), Spare parts, high-vacuum components, other types of Products replaced within the warranty period shall be handed over to Siemens to be further returned to the manufacturer no later than 14 (fourteen) days as of the date of signing the documents (certificate/order form) provided by Siemens, confirming that the respective defect was eliminated. Title to such spare part shall be transferred to Siemens.

16.Complaints

The complaint may be submitted only for those issues which have not been the subject of the device acceptance procedure performed in accordance with the terms of the contract.

In case the system does not operate properly or does not perform commands as described in operating manual contact manufacturer representative service center in RF.

Siemens Healthcare Limited Liability Company (Siemens Healthcare LLC)
115093, Moscow, Dubininskaya str. 96
Phone: (495) 737-12-52;
www.siemens.ru/healthcare

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