
Medtronic

Medtronic CoreValve LLC
1851 E. Deere Ave, Santa Ana
CA 92705, USA

Declaration Letter

Medtronic CoreValve LLC., a corporation organized under the laws of the United States, having its principal place of business at 1851 E. Deere Ave, Santa Ana, CA 92705, USA, hereby swear (or affirm) that the attached documentation:

- 1- Appendix to IFU contains 59 pages;
Are true and unaltered copies of the original document.

Name : Andrew Priest
Title: Assoc. Regulatory Affairs Specialist

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
end of the radiopaque paddle attachment, it is at the point of no recapture. The deployment knob also features a trigger, which can be engaged to make macro adjustments to the capsule position. A blue hand rest connects to the deployment knob. The end of the handle features a tip-retrieval mechanism, which can be used to withdraw the catheter tip to meet the capsule after the device has been fully deployed. The catheter packaging contains an integrated loading bath and a removable tray with 3 rinsing bowls for loading and rinsing the bioprosthesis. The integrated loading bath features a mirror, which aids in accurate placement of the bioprosthesis frame paddles during loading. In addition to these features, the device packaging is swiveled and secured to facilitate the bioprosthesis loading procedure.

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The EnVeo R Delivery Catheter System is similar to the CoreValve AccuTrak Delivery Catheter System (RC 2012/12346 from 23.08.2012) in both design and function in that it incorporates a protective deployment sheath which houses and deploys the prosthesis. Similar to the AccuTrak

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Delivery Catheter System, the EnVeo R Delivery Catheter System also features an integrated handle designed to provide the user with accurate and controlled deployment. After the EnVeo R Delivery Catheter System is placed in the vicinity of the aortic annulus, the user retracts the deployment sheath, thereby precisely deploying the valve to the desired location. In use, the deployment sheath can be partially pulled back to evaluate the valve location/function prior to fully releasing the valve. In this way, the user can make slight adjustments to the valve location if needed prior to release.

With the EnVeo R Delivery Catheter System design, the resheath and/or complete recapture feature has been developed to allow for full recapture and redeployment of the Evolut R valve. The valve can be recaptured and redeployed a maximum of two times and can be recaptured a third time to remove the valve from the patient.

In use, the outer sheath can be partially pulled back to evaluate the valve position within the annulus prior to fully releasing the valve. If the radiopaque capsule marker band has not yet reached the distal end of the paddle attachment (the user is alerted prior to this point being reached via tactile feedback from the handle), the valve can be fully resheathed/recaptured and redeployed. The flare section of the capsule expands to allow the valve to be resheathed/ recaptured, and the flex section of the capsule provides column strength to support the capsule during restheating and/or recapture while maintaining flexibility. After the valve has been fully deployed, the delivery system is retracted from the patient and discarded after use.

EnVeo R Capsule

Close up of EnVeo R Paddle Attachment within Capsule

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EnVeo R Delivery Catheter System is supplied as a single use device that is sterilized using a validated Ethylene Oxide (EtO) sterilization process.

EnVeo R Delivery Catheter System (DCS) Description  Changes

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

below shows the AccuTrak Delivery Catheter System (top) compared to the EnVeo R Delivery Catheter System (bottom). As shown, the primary functional components (i.e., guidewire lumen, micro/macro controls, flush port, stability member, 18Fr sheath capsule and inner assembly distal tip) of the delivery system remain the same between the previous generation and the proposed delivery systems. The proposed EnVeo R Delivery Catheter System includes enhancements to improve the user's experience as well as the additional features of resheathing/recapture and the EnVeo R InLine sheath.

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Comparison of AccuTrak Delivery Catheter System and EnVeo R Delivery Catheter System

Characteristic	AccuTrak Delivery Catheter System	EnVeo R Delivery Catheter System	Discussion of Change
Inner Shaft Assembly			
Inner Shaft	AccuTrak DCS has a single layer of material that makes up the inner shaft	EnVeo R DCS uses a braided inner shaft. A stainless steel braid with a single axial wire is encapsulated between a lubricious inner liner and a polyimide outer jacket	The construction of the inner shaft was changed to reduce its outer diameter to retain 18Fr overall crossing profile for the system. This reduction of OD was implemented while preserving the primary function of providing a lumen for the guidewire and isolating the tissue valve from the guidewire.
Middle Member	AccuTrak DCS has a single layer of material that makes up the middle member	EnVeo R DCS uses a Middle Member, two axially opposed stainless steel axial wires covered with stainless steel braid and laminated together by jackets.	The EnVeo R DCS middle member has been reinforced to withstand the forces due to deployment and recapture of the TAV. The axial wires run the full length of the outer shaft. The addition of axial wires provides greater columnar stiffness which also gives a more stable deployment as the capsule movement is directly linked to the handle movement as the wires do not compress or extend during loading. The updated middle member provides the tensile and compressive stiffness required to

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	layer, a Pebax tie layer, a stainless steel braided layer and the outer Nylon 12 jacket	inner liner, a stainless steel braided layer, and a Pebax outer jacket	The function of the stability member to isolate the outer member from the introducer sheath and anatomy to reduce friction and improve system accuracy remains unchanged.
InLine Sheath			
InLine Sheath	N/A-AccuTrak DCS does not have an InLine sheath	The EnVeo R DCS InLine sheath provides an option for physicians to use the system without a standard Introducer sheath	The EnVeo R DCS InLine sheath allows for the effective introduction profile to be reduced by approximately 4Fr. Physicians can use the InLine sheath instead of a standard 18Fr introducer sheath, or use in combination with a standard 18Fr introducer sheath.

EnVeo R Loading System (LS)

The LS is a system of reduction cones and tubing designed to gradually reduce the diameter of the TAV radially to an optimal diameter to facilitate manual loading of the TAV into the deployment sheath capsule of the DCS. This manual loading is performed as part of procedure preparation immediately before implantation. The loading procedure is performed while immersing the LS, TAV and distal end of the DCS in cold sterile saline.

The LS comprises the following elements:

- Inflow cone
- Tip guide tube
- Outflow cone
- Sheath guide tube

Two LS models are available, one for use with the Evolut 23R TAV and one for use with the Evolut 26R and 29R. Both LS models are identical in design and color.

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Characteristic	AccuTrak Delivery Catheter System	EnVeo R Delivery Catheter System	Discussion of Change
Inner Shaft Assembly			
Tip	The AccuTrak DCS atraumatic tip is comprised of Pebax/BaSO ₄	The EnVeo R DCS atraumatic tip is comprised of Pellethane/BaSO ₄	The tip was updated to be compatible with the overall EnVeo R system, including capsule and inner member. There is no change to the principle of operation as the tip function is identical to the AccuTrak DCS.
Inner Shaft	AccuTrak DCS has a single layer of material (Vitrex PEEK 381G) that makes up the inner shaft	EnVeo R DCS uses a braided inner shaft. A stainless steel braid with a single axial wire is encapsulated between a lubricious PTFE-polyimide inner liner and a polyimide outer jacket	The construction of the inner shaft was changed to reduce its outer diameter to retain 18Fr overall crossing profile for the system. This reduction of OD was implemented while preserving the primary function of providing a lumen for the guidewire and isolating the tissue valve from the guidewire.
Middle Member	AccuTrak DCS has a single layer	EnVeo R DCS uses a Middle	The EnVeo R DCS middle member has

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The components that required upsizing between the EnVeo R Delivery Catheter System, outer diameter 18 Fr (6 mm) design and EnVeo R Delivery Catheter System, outer diameter 20 Fr (6.7 mm) presented in this supplement are as follows:

Nitinol Capsule – The nitinol capsule is identical to the commercial device with the exception that all the subcomponents were upsized to accommodate the larger size TAV, including the capsule inner and outer diameters, which increased by 2 Fr to accommodate the larger sized valve. In addition to the overall upsizing of the subcomponents, a minor length change increased the leg length by 0.6 mm due to the angle of the capsule ribs.

Distal Tip – The Distal Tip inside the Middle Member is identical to the commercial device with the exception of the outer diameter, which was increased to fit the new 20 Fr capsule.

Spindle – The Spindle inside the Middle Member is identical to the commercial device with the exception that the outer diameter was increased to maintain the same clearance with the 20 Fr capsule that is seen with the 18 Fr system.

Outflow support – The Outflow Support inside the Middle Member is identical to the commercial device with the exception that the outer diameter was increased in line with the Spindle.

InLine Sheath (ILS) – The ILS is required to maintain the same outer profile as the Capsule, so the profile of the ILS was increased to 20 Fr in alignment with the new 20 Fr capsule. The shaft changed to a thicker wall to achieve the 20 Fr outer diameter (OD). In addition, a softer grade of Pebax® material is used in the jacket to maintain flexibility. Finally, the hub design is identical to the commercial device with the exception of the increased inner diameter (ID) to mate with the 20 Fr outer shaft.

Stylet – The dimensions of the Stylet were increased to prevent use with the smaller, currently-approved Loading System. In addition, the color of the Stylet was changed from clear to blue to further differentiate the EnVeo R 20 Fr DCS from the smaller commercial DCS.

Although the DCS profile was increased from 18 Fr to 20 Fr for the Evolut 34R TAV size, the delivery system profile to Min Vessel Diameter Ratio of 1.2 was held constant, as shown below.

- 18 Fr: 6 mm Profile / 5 mm Min Vessel = 1.2 Device-to-Vessel Diameter Ratio

- 20 Fr: 6.67 mm Profile / 5.5 mm Min Vessel = 1.2 Device-to-Vessel Diameter Ratio

Maintaining this level of device interaction with the access vessel

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anatomy is expected to deliver comparable vascular access complication rates to those of the commercially-approved 18 Fr EnVeo R DCS. EnVeo R 20 Fr DCS component materials are the same as those used in the commercial EnVeo R 18 Fr DCS with the exception of several minor differences

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Devices supplied in packaging marked with the following information:

- device name and reference number;
- serial number of the product;
- trademark and address of the manufacturer;
- Single-Use
- Sterile, EO
- CE Mark (label complies with European Council 93/42 EEC);
- date of production and expiry date;
- storage conditions;
- the characteristics of the product;
- warning and informing signs.

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- name of medical device in Russian language;
- name of the manufacturer;
- country of origin;
- code UPN (unique product number);
- code CFN (customer facing number);
- date and number of Registration Certificate;
- sign of conformity in Russia ;
- warning and informing symbols and inscriptions.

Sample of the sticker

	Product name	
	Name of the manufacturer	Code UPN
	Country of origin	Code CFN
	See IFU, Date and Number of RC	

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Medtronic
EnVeo™ R Delivery Catheter System

en Delivery Catheter System / bg Система за доставка катетер / ca Sistema de entrega de catèters / da Indføringskatersystem / de Katheterschleusensystem / el Σύστημα καθετήρα παράδοσης / es Sistema de catéter de liberación / et Saadusustsüsteem / fi Latausjärjestelmä / fr Système de cathéter de livraison / he סיסטמה של קתטר לנתיב / hu Betétrendszer / it Sistema di caricamento / it Sistema de chargement / ja 送達用カテーテルシステム / ko 공급용 카테터 시스템 / lt Įleidimo kateterio sistema / lv Ievades sistēma / mk Система за вносување / nl Laadsysteem / no Ladesystem / pl System ładowania / pt-br Sistema de carregamento / pt-pt Sistema de carregamento / ro Sistem de încărcare / ru Приспособление для загрузки катетера / sk Zavádzací systém / sl Sistem za vstavljanje / sr Sistem za ubacivanje / sv Inmatningsystem / tr Yükleme Sistemi / uk Система завантаження

7.6 cm
 107 cm
 18 Fr (6 mm)

REF **ENVEOR-L**

LOT

0.889 mm (0.035 in)

STERILE EO

CE 0066

Medtronic CardioVascular LLC
 1651 E. Deere Avenue
 Santa Ana, CA 92705
 USA
 www.medtronic.com

Medtronic Ireland
 Parkview Business Park West
 Galway, Ireland
 Made in Ireland

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

EnVeo R DCS Packaging

The EnVeo R Delivery Catheter System (DCS) is available in one configuration that accommodates all three transcatheter aortic valves sizes.

The packaging configuration for the EnVeo R DCS includes primary packaging (sterile barrier/tray), secondary packaging (carton) and a tertiary package (shipping/distribution container). The EnVeo R DCS is loaded in a thermoformed tray which is assembled and inserted in to a single Film/Tyvek pouch which is heat sealed. The labeled pouch is contained within a labeled shelf carton. Shelf cartons will be packaged in 4 pack corrugated fiberboard shipper.


EnVeo R DCS Packaging Materials Summary

Packaging Component	Material
Primary Package (Sterile Barrier System)	
Packaging Tray - proximal	0.040" White Tritan w/ ID (internal denest)
Packaging Tray - distal	0.040" White Tritan w/ ID (internal denest)
Bath Insert	0.040" White PETG w/ ID (internal denest)
Tray Clips	0.040" Blue PETG w/ ID (internal denest), Color Pantone - 293 C

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
Pouch	<i>Top Web:</i> 20 µm Orientated Polyamide (Nylon) Film /Adhesive / 20 µm Orientated Polyamide (Nylon) Film /Adhesive / 50 µm Polyethylene <i>Bottom Web:</i> <i>Tyvek 1073B Uncoated Ink: PMS Black</i>
Mirror Label	FLEXmark® MM 500 Silver PTP with V-29 adhesive
Pouch label	<i>Paper Type:</i> Fastransfer Superior S2000n /Fastransfer <i>Adhesive:</i> S2000n/S2045
Secondary Package	
Shelf Carton	E fluted (white--white) single wall corrugated fiberboard
Shelf Carton Label	<i>Paper Type:</i> Fastransfer Superior S2000n /Fastransfer <i>Adhesive:</i> S2000n/S2045
Shelf Carton Closure Label	Clear polystyrene with V88 permanent pressure sensitive rubber based adhesive
Carton Closure Label	Acetate label
Insert	INSERT, EVOLUT R, DCS, LS,
Tertiary Package (Shipping/Distribution Container)	
Shipping Case	Single wall corrugated fiberboard (burst ≥ 275 psi or ECT ≥ 44 psi)

Packaging Description EnVeo R loading system

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Packaging Component	Material	Function
Primary Package (Sterile Barrier System)		
Inner Pouch	Top Web: 20 μ m Orientated Polyamide (Nylon) Film / Adhesive / 20 μ m Orientated Polyamide (Nylon) Film / Adhesive / 50 μ m Polyethylene Bottom Web: Tyvek 1073B Uncoated Ink: PMS Black	Sterile barrier system. Contains the LS product components.
Outer Pouch	Top Web: 20 μ m Orientated Polyamide (Nylon) Film / Adhesive / 20 μ m Orientated Polyamide (Nylon) Film / Adhesive / 50 μ m Polyethylene Bottom Web: Tyvek 1073B Uncoated Ink: PMS Black	Sterile barrier system. Provides surface for labeling.
Pouch label	Top Web: 20 μ m Orientated Polyamide (Nylon) Film / Adhesive	Printed with product specific information and

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	/ 20 m Orientated Polyamide (Nylon) Film / Adhesive / 50 m Polyethylene	adhered to the pouch.
Secondary Package		
Shelf Carton	Solid bleached sulfate	Contains the pouched product and is a surface for labeling.
Shelf Carton Label	Paper label	Printed with product specific information and adhered to the pouch and carton
Shelf Carton Closure Label	Clear polystyrene with V88 permanent pressure sensitive rubber based adhesive	Provides protection to ensure that the carton remains closed during transit.

Packaging validation of the TAV, DCS, and LS demonstrated that the packaging of the three components of the Evolut R system maintained their integrity and preserved the sterility of the packaged product during simulated handling, distribution and storage.

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EnVeo R Delivery Catheter System (Catheter)


The catheter facilitates the placement of the bioprosthesis within the annulus of the aortic valve. The catheter assembly is flexible and compatible with a 0.035-in (0.889-mm) guidewire. The distal (deployment) end of the system features an atraumatic, radiopaque catheter tip and a capsule that covers and maintains the bioprosthesis in a crimped position. The capsule includes a distal flare to enable the bioprosthesis to be partially or fully recaptured after partial deployment. A stability layer is fixed at the handle and extends down the outside of the catheter shaft. It provides a barrier between the retractable catheter and the introducer sheath and vessel walls, thus enabling the catheter to retract freely. An EnVeo R InLine sheath is assembled over the stability layer, which functions as a hemostatic introducer sheath and minimizes the access site size to the capsule diameter. The EnVeo R InLine sheath is also compatible with an 18-Fr (6-mm) introducer. The delivery catheter system consists of a catheter with an integrated handle to provide the user with accurate and controlled deployment. The handle is on the proximal end of the catheter and is used to load, deploy, recapture, and reposition the bioprosthesis. The handle features a gray front grip used to stabilize the system. The deployment knob turns to deploy the bioprosthesis precisely. Arrows on the deployment knob indicate the direction of rotation required to deploy the bioprosthesis. If desired, the deployment knob can be turned in the opposite direction to partially or fully recapture the bioprosthesis if the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture. The deployment knob also features a trigger, which can be engaged to make macro adjustments to the capsule position. A blue hand rest connects to the deployment knob. The end of the handle features a tip-retrieval mechanism, which can be used to withdraw the catheter tip to meet the capsule after the device has been fully deployed. The catheter packaging contains an integrated loading bath and a

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removable tray with 3 rinsing bowls for loading and rinsing the bioprosthesis. The integrated loading bath features a mirror, which aids in accurate placement of the bioprosthesis frame paddles during loading. In addition to these features, the device packaging is swiveled and secured to facilitate the bioprosthesis loading procedure.

DCS Subassemblies


The DCS comprises the following 4 main functional subassemblies, which are described in detail in the following sections:

Handle

Inner Shaft/Distal Tip Subassembly

18 Fr Sheath Capsule Subassembly (deployment sheath)

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
Number	Component
1	Handle
2	Inner Shaft
3	12 Fr Braided Shaft (Deployment Sheath)
4	Stability Layer

Handle

The handle provides the user the ability to (1) advance the deployment sheath during the loading of the TAV into the capsule and (2) retract the deployment sheath in order to precisely deploy the TAV to the desired location. The handle has a macro control (cursor) and a micro control (thumbwheel). The macro control is used to make gross movements while the micro control is used to make fine movements to the deployment sheath. The deployment sheath can be advanced and deployed by using either of the two controls on the handle. The macro control is slid backward to retract the deployment sheath and slid forward to advance the deployment sheath forward. The micro control is rotated to make fine movements to advance or retract the deployment sheath and is primarily used to deploy the TAV.

The proximal luer connection on the handle is used to flush the guidewire lumen with saline while preparing the DCS before the implantation procedure.

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
Inner Shaft/Distal Shaft Subassembly


and a stationary attachment point for the TAV.

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Number	Component
1	Plunger Assembly
2	PEEK Inner Tube
3	0.080" PEEK Tube
4	Hypo Tube
5	Loctite Adhesive, as required

18 Fr Sheath Capsule Subassembly


The 18 Fr sheath capsule subassembly consists of the 18 Fr capsule that houses the compressed TAV and a proximal 12 Fr braided shaft that connects the 18 Fr sheath to the moving parts of the handle. The capsule has a distal radiopaque platinum/ iridium marker to allow the user to visualize the distal edge of the capsule during deployment of the TAV. The capsule is transparent in order for the user to see the compressed TAV to ensure correct loading. A purge tube is provided at the proximal end of the sheath to allow the user to purge air that may be entrapped in the capsule area around the TAV.

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
Number	Component
1	18Fr Sheath Assembly
2	Capsule
3	12Fr Braided Shaft
4	Purge Tube with Luer
5	Mounting Boss with Side Port
6	Loctite Adhesive, as required
7	Epoxy Adhesive, as required

Stability Layer

The stability member is comprised of a 15 Fr braided sheath that fits over the proximal 12 Fr braided shaft and covers approximately 80% of the length of the deployment sheath. The proximal 2 cm of the stability shaft has a thicker wall (referred to as an internal strain relief) to increase the integrity of the shaft where it transitions into the handle. An additional component (stability shaft cap) has been added to facilitate attachment of the stability shaft to the handle. Within the stability shaft cap is an ePTFE seal and seal retainer to provide hemostasis at the stability shaft cap/outer shaft interface. An additional flush port is also provided, which allows the user to purge air that may be entrapped between the stability member and the outer shaft.

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
Number	Component
1	Stability Shaft
2	Stability Member Cap
3	ePTFE Seal
4	Seal Retainer
5	Purge Tube with Luer
6	Loctite, as required

Mass EnVeO R Delivery Catheter System, outer diameter 18 Fr (6 mm); -259,5 g

Mass EnVeO R Delivery Catheter System, outer diameter 20 Fr (6.7 mm) – 260.0 g;

The components that required upsizing between the commercially-approved design and the EnVeO R 20 Fr DCS design presented in this supplement are as follows:

Nitinol Capsule – The nitinol capsule is identical to the commercial device with the exception that all the subcomponents were upsized to accommodate the larger size TAV, including the capsule inner and outer diameters, which increased by 2 Fr to accommodate the larger sized valve. In addition to the overall upsizing of the subcomponents, a minor length change increased the leg length by 0.6 mm due to the angle of the capsule ribs. An image of the 20 Fr capsule is provided in Figure below.

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


Figure: EnVeo R 20 Fr DCS Capsule

Distal Tip – The Distal Tip inside the Middle Member is identical to the commercial device with the exception of the outer diameter, which was increased to fit the new 20 Fr capsule. Refer to Figure below for an image of the Distal Tip.


Figure: EnVeo R 20 Fr DCS Distal Tip

Spindle – The Spindle inside the Middle Member is identical to the commercial device with the exception that the outer diameter was increased to maintain the same clearance with the 20 Fr capsule that is seen with the 18 Fr system. Refer to Figure below for an image of the Spindle.

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
Figure: EnVeo R 20 Fr DCS Distal Spindle

Outflow support – The Outflow Support inside the Middle Member is identical to the commercial device with the exception that the outer diameter was increased in line with the Spindle. Refer to Figure below for an image of the Outflow Support.


Figure: EnVeo R 20 Fr DCS Distal Outflow Support

InLine Sheath (ILS) – The ILS is required to maintain the same outer profile as the Capsule, so the profile of the ILS was increased to 20 Fr in alignment with the new 20 Fr capsule. The shaft changed to a thicker wall to achieve the 20 Fr outer diameter (OD). In addition, a softer grade of Pebax® material is used in the jacket to maintain flexibility. Finally, the hub design is identical to the commercial device with the exception of the increased inner diameter (ID) to mate with the 20 Fr outer shaft. Refer to Figure below for an image of the ILS.

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

Figure 0-1: EnVeo R 20 Fr DCS Inline Sheath

Stylet – The dimensions of the Stylet were increased to prevent use with the smaller, currently-approved Loading System. In addition, the color of the Stylet was changed from clear to blue to further differentiate the EnVeo R 20 Fr DCS from the smaller commercial DCS, as shown in Figure below. The blue colorant is identical to that used in the stability shaft of the commercial device.




Figure: EnVeo R 20 Fr DCS Stylet

EnVeo R Loading System (LS)

The LS is a system of reduction cones and tubing designed to gradually reduce the diameter of the TAV radially to an optimal diameter to facilitate manual loading of the TAV into the deployment sheath capsule of the DCS. This manual loading is performed as part of procedure preparation immediately before implantation. The loading procedure is performed while immersing the LS, TAV and distal end of the DCS in

		cold sterile saline. The LS comprises the following elements: • Inflow cone • Tip guide tube • Outflow cone • Sheath guide tube  The diagram shows a horizontal assembly of components. From left to right: a 'Sheath guide tube' (a long thin tube), an 'Outflow cone' (a small cone-shaped part), a 'Back plate' (a small rectangular plate), an 'Inflow' (a vertical tube with a downward arrow), and a 'Tip guide tube' (another long thin tube). Arrows point from the text labels to their respective parts in the diagram.
	EnVeo R LS 26\29cm Drawing/ EnVeo R Приспособление для загрузки (26\29см) размеры:	

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
Medical device is sterilized via Ethylene Oxide (EtO) in accordance with Medtronic internal quality control procedures, EN/ISO 11135:2007 *Medical Device – Validation and Routine Control of Ethylene Oxide Sterilization* and residual testing was conducted per ISO 10993-7:2008 *Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals*. Both pre- and post-sterile Limulus Amoebocyte Lysate (LAL) testing was conducted utilizing the requirements from current U.S Pharmacopeial USP Standard Section 85 (Bacterial Endotoxins Test). Medical device is sterilized to a sterility assurance level of 10⁻⁶. Functional testing of the EnVeo R DCS has validated a maximum of two sterilization cycles.

The medical device is intended for single use and is not intended for reuse or re-sterilization by the user.

The EnVeo R 20 Fr DCS and EnVeo R 20 Fr LS is sterilized via Ethylene Oxide (EtO) in accordance with Medtronic internal quality control procedures, EN ISO 11135-1:2007 *Medical Device – Validation and Routine Control of Ethylene Oxide Sterilization* and residual testing was conducted per EN ISO 10993-7:2008/AC:2009 *Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals*. The sterilization process is identical to the sterilization method used for the EnVeo R 18 Fr DCS and EnVeo R 18 Fr LS.

The testing was successfully completed for the EnVeo R 20 Fr DCS and EnVeo R 20 Fr LS. The testing confirmed the manufacturing process and material changes did not impact density, EtO Residual testing,

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
Declaration Letter

Medtronic CoreValve LLC, a corporation organized under the laws of the United States, having its principal place of business at 1851 E. Deere Ave, Santa Ana, CA 92705, USA, hereby swear (or affirm) that the attached documentation:


- 1- IFU contains 29 languages and 314 pages;
- Are true and unaltered copies of the original document.

Andrew Priest
Assoc. Regulatory Affairs Specialist
Medtronic CoreValve LLC

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Maria Brittany Cook 24 Oct 2017
 MARIA BRITTANY COOK
 Notary Public
 Minnesota
 My Commission Expires
 Jan 31, 2021

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Упатства за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucțiuni de utilizare • инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з експлуатації

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-
- This diagram illustrates the exploded view of the Knauf ceiling grid assembly. The components are numbered as follows:
- 1:** Main grid section with a central rectangular cutout.
 - 2:** Small circular component, likely a fastener or clip.
 - 3:** Multiple small circular components, likely fasteners or clips, used for securing the grid.
 - 4:** A rectangular panel or insert that fits into the central cutout of the main grid.
 - 5:** A long, narrow rectangular component, likely a trim or edge piece.
 - 6:** A small circular component, likely a fastener or clip.
 - 7:** A small L-shaped component, likely a fastener or clip.
 - 8:** A long, narrow rectangular component, likely a trim or edge piece, shown in a separate exploded view.

10


Figure 4. EnVeO™ R delivery catheter system proximal tray / Фигура 4. Проксимална тава на системата на доставящ катетър EnVeO™ R / Obrázek 4. Proximální podnos systému zaváděcího katétru EnVeO™ R / Figur 4. EnVeO™ R-indförsingskatersystemets proksimale bakke / Abbildung 4. EnVeO™ R Kathetereinföhrsystem, proximale Schale / Εικόνα 4. Εγγύς δίσκος συστήματος καθετήρα παράδοσης EnVeO™ R / Figura 4. Bandedja proximal del sistema de catéter de liberación EnVeO™ R / Joonis 4. Sisestuskateetri süsteemi EnVeO™ R proksimaalne alus / Kuva 4. EnVeO™ R -sisäänvientikatetrijärjestelmän proksimaalinen tarjotin / Figure 4. Plateau proximal du système de pose EnVeO™ R / Slika 4. Proksimalna plitica kateterskog sustava za uvođenje EnVeO™ R / 4. ábra. EnVeO™ R bevezetőkátéter-rendszer proximális tálcája / Figura 4. Vassoio prossimale del sistema del catetere di rilascio EnVeO™ R / 4 paveikslēlis. „EnVeO™ R” įvedimo kateterio sistemos proksimalinis dėklas / 4. attēls. EnVeO™ R katetra ievadsistēmas proksimālā paplāte / Слика 4. Проксимален послужавник на катетерскиот систем за вградување EnVeO™ R / Afbeelding 4. Proximale tray EnVeO™ R-plaatsingskathetersysteem / Figur 4. Proksimalt brett for EnVeO™ R-kateterinnförsingssystem / Rycina 4. Taca proksymalna systemu cewnika wprowadzającego EnVeO™ R / Figura 4. Bandedja proximal do sistema de cateter introdutor EnVeO™ R / Figura 4. Tabuleiro proximal do sistema de cateter introdutor EnVeO™ R / Figura 4. Tavă proximală pentru sistem de poziționare cu cateter EnVeO™ R / Рисунок 4. Проксимальный лоток системы катетерной доставки EnVeO™ R / Obrázok 4. Proximálny zásobník systému aplikácie katétra EnVeO™ R / Slika 4. Proksimalni pladenj katetrskega sistema za dovajanje EnVeO™ R / Slika 4. Proksimalno ležište EnVeO™ R sistema katetera za plasiranje / Figur 4. Proximal bricka för EnVeO™ R införsingskatersystem / Şekil 4. EnVeO™ R taşıma kateteri sistemi proksimal tablası / Рисунок 4. Проксимальный лоток катетерной системы доставки EnVeO™ R

1. Tray connector / Конектор за тава / Spojka podnosu / Bakkekonnektor / Schalen-Verbinder / Σύνδεσμος δίσκου / Conector de la bandedja / Aluse konnektor / Tarjottimen liitin / Connecteur du plateau / Priključak plitice / Tálcacsatlakozó / Connettore del vassoio / Deklo jungtis / Paplāšu savienotājs / Конектор на послужавникот / Trayconnector / Brettkobling / Łącznik tacy / Conector da bandedja / Conector do tabuleiro / Conector tavă / Конектор лотка / Konektor zásobníka / Konektor pladnja / Konektor ležišta / Brickanslutning / Tabla konektörü / Конектор лотка
2. Swivel hinge / Шарнирна панта / Otočný závěs / Kuglehængsel / Schwenkgelenk / Στροφικός Bisagra giratoria / Pöördhing / Kääntönivel / Charnière pivotante / Okretna šarka / Forgócsukló / Cardine girevole / Sukimosi ašis / Grožama vira / Шарки за въртене / Scharnier / Svinghengsel / Zawias / Dobradiça giratória / Dobradiça de rotação / Balama tip pivot / Шарнирная петля / Otočný spoj / Vrtljiv tečaj / Šarke za okretanje / Vridbar led / Döner menteşe / Шарнірна петля
3. Clip holder / Държач за скоба / Držák svorky / Klemmeholder / Clip-Halter / Συγκρατητήρας κλίπ / Soporte para clip / Klambrihoidik / Klipsin pidike / Porte-pince / Držáč kvačice / Kapocsrögzítő / Supporto della clip / Spaustuko laikiklis / Skavu turētājs / Штипка држач / Cliphouder / Klemmeholder / Uchwytyt klipsa / Suporte do clipe / Suporte do agraf / Dispozitiv de prindere a inelului / Держатель зажима / Držiak svorky / Držalo sponke / Držáč hvataljke / Klämhållare / Klips tutucu / Тримач затискача
4. Mirror / Орледало / Zrcadlo / Spejl / Spiegel / Κάτοπτρο / Espejo / Peegel / Peili / Miroir / Zrcalo / Tükör / Specchio / Veidrodis / Spogulis / Орледало / Spiegel / Speil / Lusterko / Espelho / Espelho / Oglindă / Зеркало / Zrkadlo / Ogledalce / Ogledalo / Spegel / Ayna / Дзеркало
5. Integrated loading bath / Вградена вана за зареждане / Integrovaná lázeň pro zasunutí / Integreret isætningskar / Integrierte Aufnahmewanne / Ενσωματωμένο λουτρό φόρτωσης / Baño de montaje integrado / Integreeritud laadimisvann / Kiinteä latausmalja / Bain de chargement intégré / Integrirana posuda za punjenje / Integrált betöltőtál / Bacinella di caricamento integrata / Integruota įstatymo vonelė / Iebūvēta ievietošanas tvertne / Вграден длабок сад за составување / Geïntegreerd laadbad / Integrert innsettingsbad / Zintegrowane naczynie na kąpiel przed załadowaniem / Banho de carregamento integrado / Banho de carregamento integrado / Baie de încărcare integrată / Встроенная ванночка для загрузки / Integrovaná zavádzacia vanička / Integrirana kopel za vstavljanje / Integrirana posuda za postavljanje / Integreret innmatningsbad / Dahili yùkleme banyosu / Вбудована ванночка для завантаження

6. Tray tab / Клипс на тавата / Výstupek podnosu / Bakkeflig / Schalenlasche / Γλωσσίδι δίσκου / Pestaña de la bandeja / Aluse sakk / Tarjottimen kieleke / Picot du plateau / Jezičac plitice / Tálcafűl / Linguetta del vassoio / Dėklo auselė / Paplātes izcilnis / Јазиче на послужавникот / Traylipje / Brettklaff / Występ tacy / Aba da bandeja / Aba do tabuleiro / Lamela tăvii / Ушко лотка / Půtko zásobníka / Pritrdišče na pladnju / Jezičak ležišta / Brickflik / Tabla tirnaği / Вушко лотка
7. Locking clip / Заключваща скоба / Zajišťovací svorka / Låseklemme / Arretierclip / Клип ασφάλισης / Clip de bloqueo / Lukustusklammer / Lukitusklipsi / Pince de verrouillage / Kvačica za zaključavanje / Zárókapocs / Clip di blocco / Fiksavimo sprastukas / Fiksēšanas skava / Штипка за заключување / Vergrendelclip / Låseklemme / Klips blokujący / Clipe de bloqueo / Agrafo de bloqueo / Inel de blocare / Фиксирующий зажим / Istiaca svorka / Zaklepna sponka / Hvataljka koja se zaključava / Låsklamma / Kilitleme klipsi / Фіксуєчий затискач
8. Rinsing bowls / Купи за изплакване / Proplachovací miský / Skylleskåle / Spülwannen / Λεκάνες έκπλυσης / Cubetas de enjuague / Loputuskauhid / Huuhtelumaljat / Bols de rinçage / Posude za ispiranje / Oblitótálak / Vaschette di risciacquo / Skalavimo indai / Skalošanas tvertnes / Садови за испирање / Spoelkommen / Skylleskåler / Miseczki do płukania / Bacias de lavagem / Bacias de lavagem / Vase de clătire / Емкости для промывки / Preplachovacie nádoby / Posode za izpiranje / Posode za ispiranje / Sköljskålar / Durulama çanakları / Емності для промивання
9. Tray tab holder / Държач за клипса на тавата / Držák výstupku podnosu / Bakkefligholder / Schalenlaschen-Halter / Συγκρατητήρας γλωσσίδιού δίσκου / Soporte para pestaña de la bandeja / Aluse saki hoidik / Tarjottimen kielekkeen pidike / Porte-picot du plateau / Držáč jezička plitice / Tálcafűl-rögzítő / Supporto della linguetta del vassoio / Dėklo auselės laikiklis / Paplātes izcilņa turētājs / Држач на јазичето на послужавникот / Houder voor traylipje / Brettklaffholder / Uchwył na występ tacy / Suporte da aba da bandeja / Suporte da aba do tabuleiro / Dispozitiv de prindere a lamelei tăvii / Держатель ушка лотка / Držiak půtky zásobníka / Držalo pritrdišča na pladnju / Držáč jezička ležišta / Brickflikshållare / Tabla tirnaği tutucu / Тримач вушка лотка


Figure 5. EnVeo™ R loading system / Фигура 5. Система за зареждане EnVeo™ R / Obrázek 5. Zasouvací systém EnVeo™ R / Figur 5. EnVeo™ R-indføringsystem / Abbildung 5. EnVeo™ R Ladesystem / Εικόνα 5. Σύστημα φόρτωσης EnVeo™ R / Figura 5. Sistema de montaje EnVeo™ R / Joonis 5. Laadimissüsteem EnVeo™ R / Kuva 5. EnVeo™ R -latausjärjestelmä / Figure 5. Système de chargement EnVeo™ R / Slika 5. Sustav za potiskivanje u kateter EnVeo™ R / 5. ábra. EnVeo™ R betöltőrendszer / Figura 5. Sistema di caricamento EnVeo™ R / 5 paveikslėlis. „EnVeo™ R“ įstatymo sistema / 5. attēls. EnVeo™ R ielādes sistēma / Слика 5. Систем за составување EnVeo™ R / Afbeelding 5. EnVeo™ R-laadsysteem / Figur 5. EnVeo™ R-ladesystem / Rycina 5. System ładowania EnVeo™ R / Figura 5. Sistema de carregamento EnVeo™ R / Figura 5. Sistema de carregamento EnVeo™ R / Figura 5. Sistem de încărcare EnVeo™ R / Рисунок 5. Приспособление для загрузки клапана EnVeo™ R / Obrázok 5. Zavádzací systém EnVeo™ R / Slika 5. Sistem za vstavljanje EnVeo™ R / Slika 5. EnVeo™ R sistem za ubacivanje / Figur 5. EnVeo™ R inmatningssystem / Şekil 5. EnVeo™ R yükleme sistemi / Рисунок 5. Система завантаження EnVeo™ R

1. Catheter tip guide tube / Тръбичка за водача на върха на катетъра / Vodicí trubička hrotu katétru / Guiderør til kateterspids / Führungsrohr für die Katheterspitze / Οδηγός σωλήνας άκρου καθετήρα / Tubo de guía de la punta del catéter / Kateetri otsa juhtimistoru / Katetrin kärjen ohjausputki / Tube-guide de l'extrémité du cathéter / Cijev vodilice vrha katetera / Katéterhegy-vezető cső / Tubo guida della punta del catetere / Kateterio antgalio kreipiamasis vamzdelis / Katetra gala virzošā caurulīte / Цевка-водилка за врвот на катетерот / Kathetertipgeleidebuis / Lederør for katetertupp / Rurka prowadząca końcówkę cewnika / Tubo-guia da ponta do cateter / Tubo-guia da ponta do cateter / Tub de ghidare al vârfului cateterului / Направляющая трубка кончика катетера / Vodiaca rūrka hrotu katētra / Vodilo cevke konice katetra / Севčica vodiča vrha katetera / Ledarrör för kateterspets / Kateter ucu kılavuz hortumu / Напрямна трубка кінчика катетера

2. Inflow cone / Конус на входния поток / Vtokový kónus / Indløbskegle / Einstrom-Kegel / Κώνος εισροής / Cono de entrada / Sissevoolukoonus / Suippo sisäänvirtauskärki / Cône de compression / Ulazna cijev / Beáramlási kónusz / Cono di ingresso / Pritekêjimo kūgis / Ieplūdes konuss / Воведна инка / Instroomkegel / Tilstrømningskonus / Stożek wlotowy / Cone de entrada de fluxo / Cone de entrada / Con de admisie / Входной конус / Prftokový kužel / Vhodni lijak / Ulazna kupa / Inflödeskon / İçeri akış konisi / Вхідний конус
3. Backplate / Задна пластина / Opěrná destička / Bagplade / Rückplatte / Πλάκα στήριξης / Placa trasera / Tugiplaat / Tukilevy / Plaque arrière / Stražnja pločica / Håtlap / Piastra di supporto / Galinė plokštė / Atpakalpiūsmas plėksne / Задна плоча / Achterplaat / Bakplate / Płyta oporowa / Placa de apoio / Placa de apoio / Placă posterioară / Задняя пластинка / Oporná doštička / Hrbtna ploščica / Zadnja ploča / Bakstykke / Arka plaka / Задня пластинка
4. Outflow cone / Конус на изходния поток / Výtokový kónus / Udløbskegle / Ausstrom-Kegel / Κώνος εκροής / Cono de salida / Väljavoolukoonus / Suippo ulosvirtauskärki / Cône à vis / Izlazna cijev / Kiáramlási kónusz / Cono di uscita / Ištekêjimo kūgis / Izplūdes konuss / Одводна инка / Uitstroomkegel / Utstrømningskonus / Stożek wylotowy / Cone de saída de fluxo / Cone de saída / Con de evacuar / Выходной конус / Odtokový kužel / Izhodni lijak / Izlazna kupa / Utflödeskon / Dışarı akış konisi / Вихідний конус
5. Capsule guide tube / Тръбичка за водача на капсулата / Vodicí trubička pouzdra / Guiderør til kapsel / Kapselführungsrohr / Οδηγός σωλήνας κάψουλας / Tubo de guía de la cápsula / Kapsli juhtimistoru / Kapselin ohjausputki / Tube-guide de la capsule / Cijev vodilice omotača / Tokvezető cső / Tubo guida della capsula / Kapsulės kreipiamasis vamzdelis / Kapsulas virzošā caurulīte / Цевка-водилка за капсулата / Capsulegeleidebuis / Lederør for kapsel / Rurka prowadząca kapsułę / Tubo-guia da cápsula / Tubo-guia da cápsula / Tub de ghidare al capsulei / Направляющая трубка капсулы / Vodiaca rúrka kapsuly / Vodilo cevke kapsule / Cevčica vodiča kapsule / Ledarrör för kapsel / Kapsül kılavuz hortumu / Напрямна трубка капсули

CoreValve™ Evolut™ R System

1. Device description

The Medtronic CoreValve™ Evolut™ R system is a recapturable transcatheter aortic valve implantation system, which includes the CoreValve™ Evolut™ R transcatheter aortic valve, the EnVeo™ R delivery catheter system, and the EnVeo™ R loading system.

1.1. CoreValve™ Evolut™ R transcatheter aortic valve (bioprosthesis)

The bioprosthesis is designed to replace the native or surgical bioprosthetic aortic heart valve without open heart surgery and without concomitant surgical removal of the failed valve. The support frame is manufactured from Nitinol, which has multilevel, self-expanding properties and is radiopaque. The bioprosthesis is manufactured by suturing valve leaflets and a skirt from porcine pericardium into a tri-leaflet configuration.

The bioprosthesis is processed with alpha-amino oleic acid (AOA™), which is a compound derived from oleic acid, a naturally occurring long-chain fatty acid. AOA™ is an antimineralization treatment shown to reduce both early and late valvular calcification. The bioprosthesis is available for a range of aortic annulus diameters (Table 1).

Table 1. Patient anatomical criteria

Bioprosthesis model	Size	Aortic annulus diameter	Aortic annulus perimeter ($\pi \times$ aortic annulus diameter)
EVOLUTR-23	23 mm	18 mm to 20 mm	56.5 mm to 62.8 mm
EVOLUTR-26	26 mm	20 mm to 23 mm	62.8 mm to 72.3 mm
EVOLUTR-29	29 mm	23 mm to 26 mm	72.3 mm to 81.7 mm
EVOLUTR-34	34 mm	26 mm to 30 mm	81.7 mm to 94.2 mm

1.2. EnVeo™ R delivery catheter system (catheter)

Refer to Figure 1 and Figure 2 for catheter components.

The catheter facilitates the placement of the bioprosthesis within the annulus of the aortic valve. The catheter assembly is flexible and compatible with a 0.035 in (0.889 mm) guidewire. The distal (deployment) end of the system features an atraumatic, radiopaque catheter tip and a capsule that covers and maintains the bioprosthesis in a crimped position. The capsule includes a distal flare to enable the bioprosthesis to be partially or fully recaptured after partial deployment. A stability layer is fixed at the handle and extends down the outside of the catheter shaft. It provides a barrier between the retractable catheter and the introducer sheath and vessel walls, thus enabling the catheter to retract freely. An EnVeo R InLine™ sheath is assembled over the stability layer, which functions as a hemostatic introducer sheath and minimizes the access site size to the capsule diameter. Catheter Model ENVEOR-L is compatible with an 18 Fr (6 mm) introducer sheath. Catheter Model ENVEOR-N is compatible with a 20 Fr (6.7 mm) introducer sheath.

The delivery catheter system consists of a catheter with an integrated handle to provide the user with accurate and controlled deployment. The handle is on the proximal end of the catheter and is used to load, deploy, recapture, and reposition the bioprosthesis. The handle features a gray front grip used to stabilize the system. The deployment knob turns to deploy the bioprosthesis precisely. Arrows on the deployment knob indicate the direction of rotation required to deploy the bioprosthesis. If desired, the deployment knob can be turned in the opposite direction to partially or fully recapture the bioprosthesis if the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture. The deployment knob also features a trigger, which can be engaged to make macro adjustments to the capsule position. A blue hand rest connects to the deployment knob. The end of the handle features a tip-retrieval mechanism, which can be used to withdraw the catheter tip to meet the capsule after the device has been fully deployed.

The catheter packaging contains an integrated loading bath and a removable tray with 3 rinsing bowls for loading and rinsing the bioprosthesis. The integrated loading bath features a mirror, which aids in accurate placement of the bioprosthesis frame paddles during loading. In addition to these features, the device packaging is swiveled and secured to facilitate the bioprosthesis loading procedure (Figure 3 and Figure 4).

1.3. EnVeo™ R loading system (LS)

The LS compresses the bioprosthesis into the catheter. Refer to Figure 5 for LS components.

Refer to Table 2 for system compatibility.

Table 2. System compatibility

Bioprosthesis model	Corresponding LS model	Corresponding catheter model
EVOLUTR-23	LS-ENVEOR-23	ENVEOR-L
EVOLUTR-26	LS-ENVEOR-2629	
EVOLUTR-29		
EVOLUTR-34	LS-ENVEOR-34	ENVEOR-N

2. Warnings and precautions

Carefully read all warnings, precautions, and instructions for use for all components of the system before use. Failure to read and follow all instructions or failure to observe all stated warnings could cause serious injury or death to the patient.

2.1. Warnings

- This procedure should only be performed where emergency aortic valve surgery can be performed promptly.
- Carefully verify that the relevant patient anatomical parameters are within the specifications set forth in Table 1.
- The heart team should consider the need for revascularization when selecting a treatment option for aortic valve replacement (AVR).
- Do not freeze the bioprosthesis. Freezing causes damage to the bioprosthesis. Check the freeze indicator located inside the labeled carton. If the product has been exposed to a freeze-thaw condition, the indicator vial will break, which causes dye to escape and stain the paper backing. Do not use the bioprosthesis if this or any other sign of freezing is observed.
- Package contents are supplied sterile. Do not use if package is damaged.
- The outside of the bioprosthesis container and the outer surfaces of the catheter pouch and the LS outer pouch are nonsterile and must not be placed in the sterile field.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- If a bioprosthesis and catheter have been removed from a patient, dispose of both the bioprosthesis and catheter; do not attempt to reuse either component. Both the bioprosthesis and catheter must be replaced with new sterile components.
- Package contents should not be opened until implantation and sizing are certain.
- Do not handle or manipulate the bioprosthesis with sharp or pointed objects.
- Use products before the Use By date.
- Store the bioprosthesis at room temperature.
- Carefully remove the catheter from the packaging to prevent it from kinking.
- Perform the **Bioprosthesis rinsing procedure** (Section 7.3) before loading the bioprosthesis into the catheter.
- After the bioprosthesis has been loaded into the catheter, maintain the distal end of the catheter with the bioprosthesis in cold (0°C to 8°C [32°F to 46°F]) sterile saline until implantation.
- Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.
- Mechanical failure of the delivery catheter system and/or accessories may result in patient complications.
- The safety and effectiveness of the bioprosthesis for aortic valve replacement has not been evaluated in patient populations presenting with blood dyscrasias as defined: leukopenia (WBC <1000 cells/mm³), thrombocytopenia (platelet count <50,000 cells/mm³), history of bleeding diathesis or coagulopathy, or hypercoagulable states.

2.2. Implantation precautions

- Before catheter insertion, remove the loading stylet.
- Before catheter insertion, perform balloon aortic valvuloplasty (BAV) of the native aortic valve to facilitate placement of the bioprosthesis.
- While the catheter is in the patient, ensure the guidewire is extending from the proximal end of the catheter. Do not remove the guidewire from the catheter while the catheter is inserted in the patient.
- There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture).
- For direct aortic access procedures, use caution when advancing and maintaining the introducer sheath position to avoid tissue and vascular injury and unacceptable access site bleeding.

2.3. Repositioning precautions

- If the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be recaptured or repositioned. During deployment, the deployment knob provides a tactile indication as a notification before the point of no recapture.
- Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment (point of no recapture), retrieval of the bioprosthesis from the patient (for example, use of the catheter) is not recommended. Retrieval after the point of no recapture may cause mechanical failure of the delivery catheter system, aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery.
- During deployment, the bioprosthesis can be advanced or withdrawn as long as annular contact has not been made. Once annular contact is made, the bioprosthesis cannot be advanced in the retrograde direction; recapture until the bioprosthesis is free from annular contact, and then reposition in the retrograde direction. If necessary, and the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be withdrawn (repositioned) in the antegrade direction. However, use caution when moving the bioprosthesis in the antegrade direction.

Caution: Use the handle of the delivery system to reposition the bioprosthesis. **Do not** use the outer catheter shaft.

- Once deployment is complete, repositioning of the bioprosthesis (for example, use of a snare and/or forceps) is not recommended. Repositioning of a deployed valve may cause aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery.

Do not attempt to retrieve or to recapture a bioprosthesis if any one of the outflow struts is protruding from the capsule. If any one of the outflow struts has deployed from the capsule, the bioprosthesis must be released from the catheter before the catheter can be withdrawn.

2.4. Postimplantation precautions

- Ensure the capsule is closed before catheter removal.
- When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.
- Exercise care when crossing an implanted bioprosthesis with the catheter or other devices.
- The safety and performance of a CoreValve™ Evolut™ R bioprosthesis implanted within a transcatheter bioprosthesis have not been demonstrated. However, in the event that a CoreValve™ Evolut™ R bioprosthesis must be implanted within a transcatheter bioprosthesis to improve valve function, valve size and patient anatomy must be considered before implantation of the CoreValve™ Evolut™ R bioprosthesis to ensure patient safety (for example, to avoid coronary obstruction).
- In the event that valve function or sealing is impaired due to excessive calcification or incomplete expansion, a postimplant balloon dilatation of the bioprosthesis may improve valve function and sealing. To ensure patient safety, valve size and patient anatomy must be considered when selecting the size of the balloon used for dilatation. The balloon size chosen for dilatation should not exceed the diameter of the native aortic annulus or, for surgical bioprosthetic valves, the manufacturer's labeled inner diameter. Refer to the specific balloon catheter manufacturer's labeling for proper instruction on the use of balloon catheter devices.

3. Potential adverse events

Adverse events that may be associated with the use of the CoreValve™ Evolut™ R system include, but may not be limited to, the following:

- Death
- Myocardial infarction, cardiac arrest, cardiogenic shock, cardiac tamponade
- Coronary occlusion, obstruction, or vessel spasm (including acute coronary closure)
- Cardiovascular injury (including rupture, perforation, tissue erosion, or dissection of vessels, ascending aorta trauma, ventricle, myocardium, or valvular structures that may require intervention)
- Emergent surgical or transcatheter intervention (for example, coronary artery bypass, heart valve replacement, valve explant, percutaneous coronary intervention [PCI], balloon valvuloplasty)
- Prosthetic valve dysfunction (regurgitation or stenosis) due to fracture; bending (out-of-round configuration) of the valve frame; underexpansion of the valve frame; calcification; pannus; leaflet wear, tear, prolapse, or retraction; poor valve coaptation; suture breaks or disruption; leaks; mal-sizing (prosthesis-patient mismatch); malposition (either too high or too low)/malplacement
- Prosthetic valve migration/embolization

- Prosthetic valve endocarditis
- Prosthetic valve thrombosis
- Delivery catheter system malfunction resulting in the need for additional re-crossing of the aortic valve and prolonged procedural time
- Delivery catheter system component migration/embolization
- Stroke (ischemic or hemorrhagic), transient ischemic attack (TIA), or other neurological deficits
- Individual organ (for example, cardiac, respiratory, renal [including acute kidney failure]) or multi-organ insufficiency or failure
- Major or minor bleeding that may require transfusion or intervention (including life-threatening or disabling bleeding)
- Vascular access-related complications (for example, dissection, perforation, pain, bleeding, hematoma, pseudoaneurysm, irreversible nerve injury, compartment syndrome, arteriovenous fistula, stenosis)
- Mitral valve regurgitation or injury
- Conduction system disturbances (for example, atrioventricular node block, left-bundle branch block, asystole), which may require a permanent pacemaker
- Infection (including septicemia)
- Hypotension or hypertension
- Hemolysis
- Peripheral ischemia
- Bowel ischemia
- Abnormal lab values (including electrolyte imbalance)
- Allergic reaction to antiplatelet agents, contrast medium, or anesthesia
- Exposure to radiation through fluoroscopy and angiography
- Permanent disability

4. Patient selection and treatment

4.1. Indications

The CoreValve™ Evolut™ R system is indicated for patients with symptomatic native aortic valve stenosis or a stenosed, insufficient, or combined surgical bioprosthetic valve failure necessitating valve replacement. The system is indicated for patients who are at high or greater risk for surgical aortic valve replacement OR are ≥75 years of age and at intermediate risk for surgical AVR (Society of Thoracic Surgeons operative risk score ≥4% or with an estimated hospital mortality ≥4% as assessed by the heart team).

4.2. Contraindications

The CoreValve™ Evolut™ R system is contraindicated for patients presenting with any of the following conditions:

- Known hypersensitivity or contraindication to aspirin, heparin (HIT/HITS) and bivalirudin, ticlopidine, clopidogrel, Nitinol (Titanium or Nickel), or sensitivity to contrast media, which cannot be adequately premedicated
- Preexisting mechanical heart valve in aortic position
- Ongoing sepsis, including active endocarditis

4.3. Access route considerations

Patients must present with transarterial access vessels with diameters that are either ≥5 mm when using Model ENVEOR-L or ≥5.5 mm when using Model ENVEOR-N, or patients must present with an ascending aortic (direct aortic) access site ≥60 mm from the basal plane. Patients must also present with anatomical dimensions as described in Table 1.

The risks of long-term anticoagulation and/or antiplatelet therapy should be considered. Implantation of the bioprosthesis is not recommended in patients with aortic root angulation (the angle between the plane of the aortic valve annulus and the horizontal plane/vertebrae) of >30° for right subclavian/axillary access or >70° for femoral and left subclavian/axillary access.

For subclavian access, use caution in patients with a preexisting patent Left Internal Mammary Artery/Right Internal Mammary Artery (LIMA/RIMA) graft. For direct aortic access, ensure the access site and trajectory are free of patent RIMA or a preexisting patent RIMA graft.

For direct aortic and subclavian access procedures, care must be exercised when using the tip-retrieval mechanism to ensure adequate clearance to avoid advancement of the catheter tip through the bioprosthesis leaflets during device closure.

For direct aortic access procedures, use a separate introducer sheath; do not use the EnVeo R InLine™ sheath. Maintain the EnVeo R InLine™ sheath at the proximal end of the catheter throughout the procedure.

For transfemoral access, use caution in patients who present with multiplanar curvature of the aorta, acute angulation of the aortic arch, an ascending aortic aneurysm, or severe calcification in the aorta and/or vasculature. If ≥ 2 of these factors are present, consider an alternative access route to prevent vascular complications.

4.4. Valve-in-valve considerations

Implanting a CoreValve™ Evolut™ R bioprosthesis in a degenerated surgical aortic bioprosthesis (valve-in-valve) should generally be avoided in patients with any of the following conditions:

- The degenerated bioprosthesis presents with a significant concomitant perivalvular leak (between prosthesis and native annulus), is not securely fixed in the native annulus, or is not structurally intact (for example, wireform frame fracture).
- The degenerated bioprosthesis presents with a partially detached leaflet that in the aortic position may obstruct a coronary ostium.

5. How supplied

5.1. Bioprosthesis

The bioprosthesis is sterilized with and stored in glutaraldehyde solution. The bioprosthesis is supplied sterile and nonpyrogenic in a jar with a screw cap. A freeze indicator is placed inside the labeled carton. If the freeze indicator has been activated, do not use the bioprosthesis.

5.2. Catheter

The catheter is packaged in a single-pouch configuration. The catheter is sterile if the package is undamaged and unopened. The catheter is sterilized by ethylene oxide.

5.3. LS

The LS is packaged in a double-pouch configuration. The LS is sterile if the packages are undamaged and unopened. The LS is sterilized by ethylene oxide.

6. Required materials (not included with the system)

- Sterile saline
- Sterile syringe, 10 mL (10 cc)
- Extra-support guidewire, 0.035 in (0.889 mm) × a minimum length of 240 cm
- Valvuloplasty balloon
- Standard procedure supplies

7. Instructions for use

7.1. Inspection prior to use and swivel tray setup

1. Carefully inspect the package before opening.
Caution: Do not use after the Use By date or if the integrity of the sterile package has been compromised (for example, damaged package).
Caution: Do not use the bioprosthesis if the freeze indicator has been activated.
Remove the product from the protective package.
3. Visually check that the product is free of defects. Do not use if any defects are noted.
4. Remove the locking clip attached to the rinsing bowls.
5. Remove the rinsing bowls from the integrated loading bath.
6. Remove the locking clips that connect the distal and proximal trays.
7. Lift the tray connector from the distal tray, and swivel the distal tray 180° counterclockwise.
8. Clip the tray tab on the distal tray to the tray tab holder on the proximal tray.
9. Fill the integrated loading bath with cold, sterile saline (0°C to 8°C [32°F to 46°F]).

7.2. Preparation of the catheter and LS

1. Attach a 10 mL syringe filled with sterile saline to the capsule flush port on the proximal end of the handle.
2. Carefully lift the distal end of the catheter to a near vertical orientation. To prevent kinking, do not bend the catheter severely.
3. Open the capsule and expose the paddle attachment.
Note: Use the deployment knob to open the capsule completely until the paddle attachment is fully exposed.
4. With the capsule held vertically, flush the capsule flush port. Verify that no catheter leakage is observed during any of the flushing steps. If leakage is observed, use a new system.

5. Submerge the capsule completely in the cold saline bath while flushing the capsule flush port. Continue flushing the capsule until it is completely submerged in the bath to prevent air from entering the catheter.

Note: After the bioprosthesis has been loaded into the capsule, the capsule flush port can no longer be flushed.

6. Secure a locking clip in the clip holder to angle the catheter tip into the integrated loading bath.
7. Place the LS components in the integrated loading bath.

7.3. Bioprosthesis rinsing procedure

1. Fill each of the 3 rinsing bowls (provided within the packaging) with approximately 500 mL of fresh, sterile saline at ambient temperature between 15°C to 25°C (59°F to 77°F).

Caution: Do not handle or manipulate the bioprosthesis with sharp or pointed objects. Use atraumatic forceps only.

2. Confirm the integrity of the primary bioprosthesis container. Open the container and remove the bioprosthesis. Let any remaining solution drain from the bioprosthesis completely.

Note: Retain the container with the original solution. It may be needed to store and return a rejected bioprosthesis.

3. Compare the serial number on the container with the serial number on the tag attached to the bioprosthesis.

Caution: If the serial numbers do not match, do not use the bioprosthesis.

4. Carefully remove the serial number tag from the bioprosthesis and retain the tag.
5. Immerse the entire bioprosthesis in a sterile rinsing bowl.
6. Gently agitate the bioprosthesis by hand for 15 seconds to remove the glutaraldehyde from the bioprosthesis.
7. Repeat steps 5 and 6 in one of the remaining rinsing bowls.
8. Leave the bioprosthesis submerged in sterile saline in the third rinsing bowl until it is ready to be loaded.

7.4. Bioprosthesis loading procedure

Perform the bioprosthesis loading procedure while the distal end of the catheter is immersed in the integrated loading bath filled with cold, sterile saline (0°C to 8°C [32°F to 46°F]). The bioprosthesis should remain immersed in saline during the loading process to minimize the introduction of air into the loaded system.

Note: Confirm the LS and catheter sizes are compatible with the bioprosthesis size (Table 2).

Caution: Rapid capsule advancement can contribute to difficulties with loading the valve. Slowly advancing the capsule helps facilitate successful loading.

1. Submerge and cool the bioprosthesis in the integrated loading bath filled with cold, sterile saline.
2. Advance the capsule guide tube over the catheter shaft toward the handle until the flexible tip is completely proximal to the paddle attachment and the end of the capsule is even with the edge of the rigid portion of the capsule guide tube.

Caution: Do not attempt to advance the flexible tip of the capsule guide tube over the capsule; this will prevent the capsule flare from expanding fully and prevent proper loading.

3. Ensure that the backplate has been inserted into the inflow cone and the exposed part of the backplate is facing up.
4. Insert the inflow portion of the bioprosthesis frame into the inflow cone. Ensure that the bioprosthesis frame paddles are aligned with the paddle attachment pockets.
5. Secure the outflow cone onto the inflow cone until it locks.
6. Insert the catheter tip guide tube completely into the distal end of the inflow cone. Inspect the outflow struts of the bioprosthesis and, if needed, manually manipulate so that they are evenly spaced and the bioprosthesis frame paddles are approximately 180° apart.
7. Insert the distal catheter tip into the catheter tip guide tube.

Note: Allow the loading tool to rest on the loading bath floor to ensure coaxial alignment with the catheter to assist in seating the bioprosthesis frame paddles within the paddle attachment pockets.

8. Retract the catheter tip guide tube to set the bioprosthesis frame paddles into the paddle attachment pockets.

Note: If the bioprosthesis frame paddles do not seat properly within the paddle attachment pockets upon retracting the catheter tip guide tube, slightly manipulate the position of the loading tool until paddle seating is achieved.

Note: If necessary, it is acceptable to manually compress the bioprosthesis frame paddles with fingertips to help seat the paddles within the paddle attachment pockets.

Note: Ensure both bioprosthesis frame paddles are completely seated within the paddle attachment pockets before continuing to the next step.

9. Hold the loading tool stationary with one hand, and with the other hand manually advance the capsule guide tube so that the flexible section covers the paddle attachment pockets and the top portion of the outflow struts. Use the mirror to ensure

that both bioprosthesis frame paddles are positioned correctly in the paddle attachment pockets and the outflow struts are within the flexible tip.

10. Advance the capsule to cover the bioprosthesis frame paddles, pausing when the capsule covers the proximal half of the paddles to confirm the paddles are both still properly seated before advancing further. Use the mirror to ensure that both paddles are captured in the capsule.

Caution: Do not advance the capsule over the bioprosthesis frame paddles unless they are fully seated in the center of the paddle attachment pockets. Advancing the capsule before the paddles are fully seated could damage the capsule and result in emboli.

11. Advance the capsule to capture the bioprosthesis outflow struts. Use the mirror to ensure that all bioprosthesis outflow struts are symmetrical and captured in the capsule.
12. Continue to advance the capsule until it reaches the distal end of the commissure pad of the bioprosthesis. The capsule should completely cover the commissure pad.
13. Remove the backplate and the catheter tip guide tube from the outflow cone.
14. While holding the capsule guide tube stationary, advance the inflow cone to crimp the inflow portion of the bioprosthesis frame until the outflow cone contacts the capsule guide tube.

Note: Ensure the bioprosthesis frame axis is visually aligned (coaxial) with the inflow cone axis during the insertion of the bioprosthesis into the inflow cone. Complete the insertion of the bioprosthesis into the inflow cone in one uninterrupted movement.

15. Advance the capsule over the bioprosthesis until the capsule comes within 5 mm of the catheter tip.

16. Remove the outflow cone and inflow cone from the catheter.

17. Advance the capsule to close the gap between the capsule and catheter tip completely.

Caution: Stop advancing the capsule once the gap to the catheter tip is closed. Advancing the capsule farther could damage the capsule.

18. Remove the capsule guide tube from the catheter. Slightly rotate the deployment knob in the direction of the arrows to relieve stress. Ensure that the capsule does not separate from the catheter tip.

Note: After the bioprosthesis has been loaded into the capsule, the capsule flush port can no longer be flushed.

19. Visually and tactilely inspect the capsule for a misloaded bioprosthesis. The capsule should be straight, smooth, and free of any bends, protrusions, or discolorations. If any of these conditions are felt or observed, the bioprosthesis is likely to be misloaded.

Note: If a misload is detected, unsheath the bioprosthesis and examine the bioprosthesis for damage (for example, permanent frame deformation, frayed sutures, or valve damage). Do not attempt to reload a damaged bioprosthesis; if no issues are found, a second attempt may be made to load an undamaged bioprosthesis. However, the catheter, LS, loading tray, and saline must be replaced with new sterile components. Do not load the bioprosthesis onto the catheter more than 2 times or after it has been inserted into a patient.

20. Attach a 10 mL syringe filled with sterile saline to the stability layer flush port on the distal end of the handle and flush.
21. Remove the loading stylet from the guidewire lumen at the capsule.

22. Attach a 10 mL syringe filled with sterile saline to the wire lumen flush port on the proximal end of the handle and flush.

23. Attach a 10 mL syringe filled with sterile saline to the EnVeo R InLine™ sheath flush port and flush.

24. Before inserting into a patient, visually inspect the loaded bioprosthesis under fluoroscopy.

Note: If a misload is detected, unsheath the bioprosthesis and examine the bioprosthesis for damage (for example, permanent frame deformation, frayed sutures, or valve damage). Do not attempt to reload a damaged bioprosthesis; if no issues are found, a second attempt may be made to load an undamaged bioprosthesis. However, the catheter, LS, loading tray, and saline must be replaced with new sterile components. Do not load the bioprosthesis onto the catheter more than 2 times or after it has been inserted into a patient.

25. Leave the bioprosthesis submerged in sterile saline until implantation.

7.5. Bioprosthesis implantation

Note: Use systemic anticoagulation during the implantation procedure based on hospital or physician protocol. If heparin is contraindicated, consider an alternative anticoagulant.

7.5.1. Bioprosthesis deployment

1. Prepare the vascular access site according to standard practice.
2. Predilate the native aortic valve with an appropriate diameter valvuloplasty balloon.

Note: Balloon predilatation of a stenotic surgical aortic bioprosthesis has not been evaluated. In cases where there is severe stenosis, predilatation of the surgical aortic bioprosthesis may be performed, and the steps used are identical to native valve predilatation.

3. Backload the catheter onto the guidewire. Insert the catheter tip and capsule through the access site, while maintaining the EnVeo R InLine™ sheath tip against the proximal end of the capsule. Then, insert the EnVeo R InLine™ sheath through the access site, maintaining contact with the capsule. Maintain guidewire position across the aortic valve.

Note: For direct aortic access procedures, use a separate introducer sheath; do not use the EnVeo R InLine™ sheath. Maintain the EnVeo R InLine™ sheath at the proximal end of the catheter throughout the procedure.

4. Under fluoroscopic guidance, advance the catheter over the guidewire to the aortic annulus. **Do not** rotate the catheter as it is advanced; rotating the handle does not rotate the capsule.

Caution: There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture).

Caution: Persistent force on the catheter can cause the catheter to kink, which could increase the risk of vascular complications (for example, vessel dissection or rupture).

5. Position the catheter so that the bioprosthesis is at the optimal depth relative to the valve annulus. For surgical bioprosthetic valves, consider the features of the valve when determining the optimal placement of the bioprosthesis.
6. To deploy the bioprosthesis, rotate the deployment knob in the direction of the arrows. The capsule retracts and exposes the bioprosthesis. Continue deploying the bioprosthesis in a controlled manner, adjusting valve position as necessary and noting the position of the radiopaque capsule marker band and paddle attachment.

Warning: Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.

Note: Slight antegrade repositioning of a partially deployed bioprosthesis (before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment) can be achieved by carefully withdrawing the catheter.

Caution: Use the catheter handle to reposition the bioprosthesis. **Do not** use the outer catheter shaft.

7. Before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, evaluate the bioprosthesis position.

Note: When the bioprosthesis is approximately 2/3 deployed, the deployment knob provides a tactile indication as a notification before the point of no recapture. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture.

8. Either complete bioprosthesis deployment or initiate bioprosthesis recapture.

7.5.2. Bioprosthesis recapture (optional)

The bioprosthesis is recapturable during deployment before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment. Deployment of the bioprosthesis can be attempted 3 times. If the bioprosthesis is recaptured a third time, it must be removed from the patient.

1. Rotate the deployment knob in the opposite direction of the arrows to recapture the bioprosthesis. A partially recaptured bioprosthesis can be repositioned or fully recaptured.

Warning: Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.

2. To fully recapture the bioprosthesis, continue rotating the deployment knob until the gap between the capsule and catheter tip is closed.

Caution: Stop advancing the capsule once the gap to the catheter tip is closed. Advancing the capsule farther could damage the capsule.

3. Reposition the recaptured bioprosthesis at the optimal depth relative to the valve annulus. For surgical bioprosthetic valves, consider the features of the valve when determining the optimal placement of the bioprosthesis.
4. Redeploy the bioprosthesis (Section 7.5.1, steps 6 and 7).
5. Either complete bioprosthesis redeployment or initiate bioprosthesis recapture. If the bioprosthesis has been recaptured 3 times, withdraw the recaptured bioprosthesis.

7.5.3. Catheter removal

1. Under fluoroscopic guidance, confirm that the catheter tip is coaxial with the inflow portion of the bioprosthesis.
2. For transfemoral access, withdraw the catheter until the catheter tip is positioned in the descending aorta. For direct aortic access and subclavian access, withdraw the catheter until the catheter tip is close to the distal tip of the introducer sheath.

3. Under fluoroscopic guidance, close the catheter capsule before withdrawal.

Caution: Close the capsule until it is aligned with the catheter tip. Do not overcapture the catheter tip, because it could interfere with catheter withdrawal through the introducer sheath or cause vessel trauma upon removal.

Caution: Ensure the capsule is closed before catheter removal.

Caution: When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.

4. Withdraw the catheter until the capsule meets the distal end of the EnVeo R InLine™ sheath.

Note: For direct aortic access procedures, maintain the EnVeo R InLine™ sheath at the proximal end of the catheter.

5. Withdraw the catheter and EnVeo R InLine™ sheath together, and dispose of the device in accordance with local regulations and hospital procedures.

6. Perform routine aortogram to assess the bioprosthesis for proper expansion and function.

Note: In the event that valve function or sealing is impaired due to excessive calcification or incomplete expansion, a postimplant balloon dilatation of the bioprosthesis may improve valve function and sealing. To ensure patient safety, valve size and patient anatomy must be considered when selecting the size of the balloon used for dilatation. The balloon size chosen for dilatation should not exceed the diameter of the native aortic annulus or, for surgical bioprosthetic valves, the manufacturer's labeled inner diameter. Refer to the specific balloon catheter manufacturer's labeling for proper instruction on the use of balloon catheter devices.

8. Training

Medtronic provides detailed operator training. For training material and in-service support, contact a Medtronic representative.

Caution: Implantation using the CoreValve™ Evolut™ R system should be performed only by physicians who have received CoreValve™ Evolut™ R training.

9. MRI safety information

Nonclinical testing and modeling have demonstrated that the CoreValve™ Evolut™ R bioprosthesis is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 2500 gauss/cm: (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the CoreValve™ Evolut™ R bioprosthesis is expected to produce a maximum in vivo temperature rise of less than 4.0°C after 15 minutes of continuous scanning. Based on nonclinical data, the image artifact caused by the device will extend no greater than 7 mm from the CoreValve™ Evolut™ R bioprosthesis when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Scanning under the conditions defined above may be performed immediately after implantation.

For deployment of a CoreValve™ Evolut™ R bioprosthesis inside of a failed surgical bioprosthetic aortic valve, consult the MRI modeling pertaining to the failed valve for additional image artifact information. The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

10. Disclaimer of warranty

The following disclaimer of warranty applies to customers outside the United States:

ALTHOUGH THE MEDTRONIC COREVALVE™ EVOLUT™ R SYSTEM, COREVALVE™ EVOLUT™ R TRANSCATHETER AORTIC VALVE, ENVEO™ R DELIVERY CATHETER SYSTEM, AND ENVEO™ R LOADING SYSTEM, HEREFTER REFERRED TO AS "PRODUCT", HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF

WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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© 2016 Medtronic
Trademarks may be registered and are the property of their respective owners.