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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SHENZHEN
COMEN MEDICAL INSTRUMENTS CO., LTD. on the
DECLARATION is genuine.

China Council for the Promotion
of International Trade

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Authorized Signature: Huang Chunxuan

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This is to certify that : this is user manual for Multi-parameter patient monitor C90.

This document is used for MOH registration in Russia.

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Manufacturer: Shenzhen Comen Medical Instruments Co., Ltd.

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Address: South of Floor 7, Block 5, 4th Industrial Area of Nanyou, Nanshan District, Shenzhen 518052, Guangdong, P.R. China.

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
2018

Copyright

Shenzhen Comen Medical Instruments Co.,Ltd Edition: A

Product Name: Patient Multi-parameter Monitor

Product Model: C90

Executive Standard No.:

Standard Registration No.

Product Registration No.

Production License No.

Registered Address: South Side of Floor 7 Bldg. 5, Floor 6 Bldg. 4, Fourth Industrial Zone of Nanyou, Nanshan Avenue, Nanshan District, Shenzhen City

Production Address: Floor 6 Bldg. 4, Fourth Industrial Zone of Nanyou, Nanshan Avenue, Nanshan District, Shenzhen City

Statement

Shenzhen Comen Medical Instruments Co., Ltd (“Comen” or “Comen Company”) possesses the copyright of the non-published Instruction Manual, and has the right to handle it as confidential data. This Instruction Manual is only used as the reference data for the operation, maintenance and servicing of Comen products.

The others have no right to disclose the content of this Instruction Manual to any third person.

This manual contains the proprietary information which is protected by copyright law. With copyright reserved, any person shall not photocopy or duplicate any part of this manual or translate it into any other language without the prior written permission of Comen Company.

Comen Company shall not assume responsibility for the incidental or consequential damages resulting from the errors in this manual or from the provision, actual performance and use of this manual. Comen Company shall not provide the proprietary right vested by patent law to any third party. Comen Company shall not assume legal responsibility for the legal consequences resulting from violating the patent law and the rights of any third party.

The content in this manual may be subject to change without further notice.

Guarantee

Under the following conditions, Comen Company shall be responsible for the safety, reliability and performance of the product:

- The product shall be used according to the Instruction Manual.

The installation, maintenance and upgrading of the product shall be conducted by the personnel recognized or authorized by Comen Company.

- The storage environment, working environment and electrical environment of the product comply with the product specifications.
- The serial number tag or manufacturing mark of the product shall be clear and identifiable. It is verified through the appraisal of Comen Company that this product is manufactured by Comen Company.
- The damages are caused by non-human factors (e.g. accidental fall and deliberate destruction).

The products that are under the warranty service regulations of Comen Company shall enjoy free service. As for the products that are beyond the warranty service regulations, Comen Company shall charge for the service. If the products are transported to Comen Company for maintenance, the user shall bear the freight (including the customs expenses).

Return of products

If the products need to be returned to Comen Company, please follow the following steps:

Acquisition of the right to return the goods: Contact the after-sales department of Comen Company, tell it the number of the instrument made by Comen; this number was printed on the equipment nameplate; if this instrument number is not legible, the goods returned shall not be accepted. Please specify the instrument number and production date, and briefly describe the reason for returning the goods.

After-sales Service Unit

Name: After-sales Service Department of Shenzhen Comen Medical Instrument Co., Ltd

Address: Floor 6 Bldg. 4, Fourth Industrial Zone of Nanyou, Nanshan Avenue, Nanshan District, Shenzhen City

Tel.: 0755—26431236

Fax: 0755—26431232

Customer Service Hotline: 4007009488

Postal Code: 518052

Preface

In this Instruction Manual, the performance, operation method and other safety information of the C90 Patient Multi-parameter Monitor (hereinafter referred to as Monitor or Patient Multi-parameter Monitor) are described in detail. Here is the best starting point for new users to use their monitors.

Performance, Structure and Composition

The monitor is mainly comprised of host machine and corresponding functional accessories (ECG cable, noninvasive blood pressure cuff, blood oxygen sensor, body temperature sensor, IBP sensor, CO2 sensor, anesthetic module analyser, noninvasive cardiac output chief cable, noninvasive cardiac output lead, invasive cardiac output float catheter, invasive cardiac output injecta temperature probe, anesthetic depth sensor, and anesthetic depth patient cable, etc.).

Applicable Scope

The Patient Multi-parameter Monitor ("the Monitor") is applicable to the monitoring or measuring of patient physiological sign parameters such as ECG, noninvasive blood pressure, pulse oxygen saturation, body temperature, respiration, pulse frequency / pulse, invasive blood pressure, carbon dioxide of the expiration tip, anesthetic gas, noninvasive cardiac output, invasive cardiac output and anesthetic depth according to the relevant parameters of module dynamic plugging automatic identification, with the functions of respiratory oxygenation chart analysis, ST segment analysis, arrhythmia analysis and drug calculation. The monitor can be used to monitor the adult, child and newborn in the operating room, ICU, CCU and bedside.

Applicable Objects of the Product

This instruction manual is applicable to the professional clinical medical staff or the persons who are experienced in using the monitoring equipment for reading. The readers shall have the knowledge and working experience in medical procedure, practice and terms necessary for monitoring the patients.

Figures

All the figures provided in this Instruction Manual shall be for your reference only. The menus, options, values and functions in the figures may be not entirely consistent with what you see from the monitor.

Conventions

- —>: This symbol is used to indicate the operating step.
- Character is used to indicate the character string in the software.

Content

Chapter 1 General	1
Safety Information	1
Overview of the Monitor	3
Screen Display	4
Key Functions and Basic Operations	7
Key Functions	7
Basic Operations	7
External Interface of the Monitor	8
Left Panel of the Monitor	8
Right Panel of the Monitor	9
Back Cover of the Instrument	10
Equipment Symbols	11
Chapter 2 Installation of the Monitor	12
2.1 Unpacking and Examination	12
2.2 Connect the AC power cord	12
2.3 Power on	12
2.4 Connecting the Sensor	13
Chapter 3 System Settings	14
3.1 Measuring Settings	14
3.2 Main Menu	15
3.2.1 Patient Management	15
3.2.2 Alarm Settings	18
3.2.3 Recording Output Settings	18

3.2.4	Review	18
3.2.5	Drug Calculation.....	19
3.2.6	Interface Switch	19
3.2.7	Monitor Information	19
3.2.8	Monitor Maintenance (user maintenance)	19
3.2.9	Demo Features	20
3.2.10	System Time Settings.....	20
3.2.11	Network Bed Number	21
Chapter 4 System Working Interface.....		22
4.1	System Working Interface.....	22
4.1.1	Normal Interface	22
4.1.2	Large Font Interface.....	23
4.1.3	Tendency Coexistence Interface	23
4.1.4	Respiratory Oxygenation Chart Interface	25
4.1.5	List Interface	26
4.1.6	Other Bed Observation	26
4.1.7	Other Interfaces.....	27
4.2	Setting the Interface Style	28
4.2.1	Setting the waveform scanning speed.....	28
4.2.2	Setting the waveform display pattern.....	28
4.2.3	Setting the display color.....	29
4.2.4	Setting the waveform to be exchanged	29
Chapter 5 Alarm		30
5.1	Alarm Type	30

5.2	Alarm Level	30
5.3	Alarm Mode	31
5.4	Lighting Alarm	31
5.5	Sound Alarm	31
5.6	Alarm	32
5.7	Alarm Parameter Flashing	32
5.8	Icon of Alarm Status	32
5.9	Setting of Alarm	33
5.10	Set the Alarm Volume	34
5.11	Setting parameter alarm	34
5.11.1	Setting the upper and lower limits of alarm as well as the ON and OFF alarms	34
5.11.2	Setting the alarm level	35
5.11.3	Alarm Record	35
5.12	Verifying the Alarm System	35
5.13	Time-out or Disenabled Alarm	36
Chapter 6 Battery		38
6.1	General	38
6.2	Installing Battery	39
6.3	Optimization and check of battery performance	39
6.4	Battery Recovery	41
Chapter 7 Data Review		42
7.1	NIBP Measurement Review	42
7.2	Warning Event Review	43
7.3	Arrhythmia Review (no window display)	44

7.4	Review of tendency chart.....	45
7.5	Tendency Table Review	47
7.6	Waveform holographic Review	48
Chapter 8 Maintenance, Cleaning and Care		50
8.1	General.....	50
8.2	Maintenance and check.....	51
8.3	Maintenance schedule.....	52
8.4	Cleaning and Sterilization.....	52
8.4.1	Note on Cleaning and Sterilization	52
8.4.2	Cleaning.....	53
8.4.3	Sterilization.....	54
Chapter 9 Patient Safety		55
9.1	Safety Instruction.....	55
9.2	Environment.....	55
9.3	Power Requirement.....	55
9.4	Grounding protection.....	55
9.5	Equipotential Grounding.....	56
9.6	Condensation	56
9.7	Description of Symbols on Instrument	57
Chapter 10 ECG Monitoring.....		58
10.1	Definition of ECG monitoring.....	58
10.2	Precautions for ECG Monitoring	59
10.3	Steps of Monitoring	60
10.3.1	Preparation	60

10.3.2	Installing ECG lead.....	61
10.4	ECG display and screen hotkey	67
10.5	Setting of ECG.....	69
10.5.1	Set the type of lead.....	69
10.5.2	Set the name of focus monitoring lead.....	70
10.5.3	Set the gain.....	70
10.5.4	Set the filtering method.....	71
10.5.5	Set the calculation channel.....	71
10.5.6	Set the power frequency inhibition	71
10.5.7	ECG calibration	72
10.6	ST segment analysis.....	72
10.6.1	About ST segment analysis.....	72
10.6.2	The impact for ST-segment.....	73
10.6.3	On/Off ST segment analysis	73
10.6.4	Point of ST segment analysis.....	73
10.7	Arrhythmia analysis.....	74
10.7.1	On/Off arrhythmia analysis.....	75
10.7.2	Arrhythmia alarm settings.....	75
10.7.3	Arrhythmia Review.....	76
Chapter 11 RESP Monitor		78
11.1	RESP (Breathing) Measurement.....	78
11.1.1	Generation of Breath.....	78
11.1.2	The Settings of Breathing Monitor	78
11.1.3	Place Electrodes for the Breathing Measurement	79

11.2	RESP Settings	80
11.2.1	Gain settings.....	80
11.2.2	Set asphyxia warning time	80
11.3	Maintenance and Cleaning.....	80
Chapter 12 SpO ₂ Monitoring		82
12.1	Definition of SpO ₂ Monitoring	82
12.1.1	Principle Of Measuring SpO ₂ Plethysmography Parameter	82
12.1.2	Recognition of the type of blood-oxygen.....	83
12.1.3	Monitoring of Blood-Oxygen/ Pulse.....	83
12.2	Attentions on Monitoring of SpO ₂ / Pulse	84
12.3	Steps of Monitoring	85
12.4	Limits in Measurement	87
12.5	SpO ₂ settings.....	88
12.5.1	Set intelligent alarm	88
12.5.2	Set average time	90
12.5.3	Set calculation sensitivity.....	90
12.5.4	Set signal indication.....	90
12.6	Maintenance and Cleaning.....	91
12.7	Masimo Information	91
Chapter 13 NIBP Monitoring		93
13.1	General.....	93
13.2	NIBP Monitoring	94
13.2.1	NIBP measurement	94
13.2.2	Action Prompt.....	96

13.2.3	Measurement Restrictions.....	97
13.2.4	NIBP Display	98
13.3	NIBP Settings.....	99
13.3.1	Setting the Measurement Mode	99
13.3.2	Setting the Interval Time.....	99
13.3.3	Consecutive Measurement	100
13.3.4	Setting the Pressure Unit.....	100
13.3.5	Setting the Pre-Inflation Value	100
13.4	NIBP Reset.....	101
13.5	NIBP Pressure Calibration	101
13.6	NIBP Air Leakage Check.....	102
13.7	Maintenance and Cleaning.....	103
Chapter 14 TEMP Monitoring.....		105
14.1	Temperature Monitoring	105
14.1.1	Temperature Measuring Setting	105
14.2	TEMP Unit Setting.....	106
14.3	TEMP Display	106
14.4	Maintenance and Cleaning.....	106
Chapter 15 CO ₂ Monitoring.....		108
15.1	General.....	108
15.2	Understanding CO ₂ Modules	109
15.3	Measuring Principle and Working Process	110
15.4	Operating Instruction for CO ₂ Connection	110
15.5	CO ₂ Measurement Procedures	112

15.6	Measuring Procedure of RESPIRONICS Branded Mainstream and Sidestream Modules.....	112
15.7	Measuring Procedure of PHASEIN Branded Sidestream and Mainstream Analyzers	115
15.7.1	Measurement Steps	116
15.7.2	Check before use.....	117
15.8	CO ₂ Setting	118
15.8.1	Turning on/off CO ₂	118
15.8.2	Setting the Pressure Unit.....	118
15.8.3	Setting Gas Compensation	118
15.8.4	Setting the Suffocation Alarming Time.....	118
15.9	Discharging Waste Gases.....	119
15.10	Maintenance and Cleaning of RESPIRONICS Branded Mainstream and Sdestream Modules	119
15.10.1	Common Cleaning	120
15.10.2	Airway Adapter for Cleaning Reusable Mainstream Sensor	120
15.10.3	Method for Sterilizing Reusable Adapter	120
15.10.4	Sterilization Times for Reusable Airway Adapter	120
15.10.5	Zeroing.....	120
15.11	PHASEIN Branded Mainstream and Sidestream Analyzer Related Information	121
15.11.1	Zeroing.....	121
15.11.2	CO ₂ Module Lighting Information	121
15.11.3	Others.....	122
Chapter 16	AG Monitoring.....	123
16.1	General.....	123
16.2	Understanding the AG Module	123
16.3	Measurement Principle of Anaesthetic Gas	124

16.4	Anesthetic Gas Display	125
16.5	MAC (Minimum Alveolar Concentration) Calculation	126
16.6	Oxygen Sensors	126
16.7	Measurement Steps and Inspection Prior to Use	127
16.8	Gas Settings	130
16.9	Setting AG.....	130
16.9.1	Setting the Anesthetic Gas Working Mode	130
16.9.2	Setting Gas Compensation	131
16.9.3	Zeroing.....	131
16.9.4	Setting the Pressure Unit.....	132
16.9.5	Setting the Suffocation Alarming Time.....	132
16.10	Anesthesia Module Lighting Information.....	133
16.11	Adverse Effects on Performance.....	133
16.12	Safety Alarm Information	134
16.12.1	ISA Sidestream Gas Analyzer Safety Warning Information	135
16.12.2	ISA Mainstream Gas Analyzer Safety Warning Information.....	136
16.13	Airway Obstruction.....	138
16.14	Discharging Waste Gases	139
16.15	Consumables.....	139
16.16	Safety Symbol Information.....	139
16.17	Patents and Trademarks	141
16.18	Maintenance.....	142
16.19	Cleaning the Analyzer.....	142
Chapter 17 IBP Monitoring		144

17.1	General.....	144
17.2	Understanding IBP Modules.....	144
17.3	Caveats for IBP Monitoring.....	145
17.4	Steps of Monitoring.....	146
17.5	Setting IBP.....	147
17.5.1	Changing the Pressure Designation	147
17.5.2	Setting the Filtering Mode	148
17.5.3	Setting Pressure Scales	148
17.5.4	Setting the Pressure Unit.....	149
17.5.5	IBP Zeroing.....	149
17.5.6	Calibration	149
17.6	Maintenance and Cleaning.....	151
17.6.1	Washing of Repeatedly Usable Pressure Sensors	151
17.6.2	Sterilization	152
Chapter 18 ICG Monitoring.....		154
18.1	General.....	154
18.2	Understanding the ICG Module.....	154
18.3	Safety Information	154
18.4	ICG Parameters.....	155
18.4.1	Hemodynamic Values	155
18.4.2	List of Hemodynamic Parameters.....	156
18.5	Noninvasive Cardiac Output Display	156
18.6	Affecting Factors	156
18.7	Steps of Monitoring	156

18.7.1	Skin Preparation.....	157
18.7.2	Installing the Sensors	157
18.8	Setting the Noninvasive Cardiac Output (ICG)	158
18.8.1	Patient Information Setting	158
18.9	Malfunctions	159
18.10	Maintenance.....	159
Chapter 19	C.O. Monitoring	161
•	Machine Profile.....	161
•	Measurement Principle	161
•	Understanding the C.O. Module	161
•	Setting the Cardiac Output (C.O.).....	162
	Setting the C.O. Coefficient.....	162
	Setting the Injecta Temperature Source	162
	Setting the Injecta Temperature Source	162
	Setting the Injecta Temperature	163
	Setting the Temperature Unit	163
	Setting the Time Interval.....	163
	Setting the Injecta Volume	163
•	Cardiac Output Measurement	164
•	Measuring the Blood Temperature.....	168
•	Affecting Factors	169
•	Hemodynamic Calculation.....	170
	Input Parameters (译者注：请不要翻译下表中的中文名称).....	170
	Input Parameters	171

- Maintenance..... 172
- Chapter 20 IoC Monitoring..... 174
 - General..... 174
 - Understanding IoC Module..... 174
 - Safety Information 174
 - IoC Display 175
 - IoC Connection 176
 - Electrode Impedance Test Results Window 178
 - IoC Setting 179
 - Set Gain..... 179
 - IoC Display Mode Setting 179
 - Trend Time Setting 179
 - Cleaning and Maintenance..... 180
 - Malfunctions 180
- Appendix I: Annex 1
- Appendix II Product Specification 7
- Appendix III System Alarm Information 25

Chapter 1 General

Safety Information

Warning

- To warn you of the conditions where serious consequence, disadvantageous matters or danger may occur. Failure to comply with the warning will result in severe personal injury or death of the user or the patient.

Caution

- To indicate potential danger or unsafe operation. If not avoided, it may lead to mild personal injury, product malfunction, damages or property loss. It may also give rise to more severe harm.

Attention

- It emphasizes primary warnings or provides descriptions or explanations so that this product can be used in a better way.

Warning

- This monitor is used for monitoring the clinical patients, so only the doctors and nurses who are qualified through training can use this monitor.
- Before use, the user shall check whether this instrument and its accessories can work normally and safely.
- The alarm volume and upper and lower limits for alarm shall be set for different patients. When a patient is monitored, the audible alarm system cannot be merely depended on. If the alarm volume is set too low or is completely turned off, the alarm will fail and the patient safety will be endangered. The most reliable patient monitoring method shall be to closely monitor the actual clinical situation of the patient.
- This instrument can only be connected to a power socket with protective grounding. If the power socket is not connected to a grounding conductor, do not use this socket, but use the rechargeable batteries for power supply.
- Do not open the shell of this instrument to avoid the possible electric shock hazard. The maintenance and upgrading of this monitor must be conducted by the service personnel trained and authorized by Comen Company.
- The disposal of packaging materials shall comply with the local laws and regulations or the waste

disposal rules and regulations of the hospital. The packaging materials must be placed away from the children.

- Do not use this instrument at the place where there are flammable articles such as anesthetic to prevent explosion or fire from happening.
- Please carefully install the power lines and the cables for various accessories to avoid the patient from being constricted or suffocated or the cables from getting entangled and keep the patient free from electrical interference.
- Do not use mobile phone near the monitor, because the mobile phone will generate a very strong radiation field and disturb the functions of the monitor.
- For the patient with pacemaker, cardiometer might measure the heart rate by the pulse of pacemaker when cardiac arrest or arrhythmia. Do not completely rely on the alarm of cardiometer. The patient with pacemaker shall be closely monitored. For the inhibiting capacity of relevant equipment on the pacemaker, refer to the Instruction Manual.
- The operators shall not touch the patients, tables and instruments during the defibrillation period.
- Before reusing these cables, check whether the function is normal.
- The equipment connected with the monitor shall form an equipotential body (the protective grounding wire is effectively connected).
- When the monitor is shared with the electrosurgery unit, the user (doctor or nurse) shall ensure the patients safety.
- The physiological waveforms, physiological parameters and alarm information, etc. displayed by this monitor shall be for the doctors' reference only and cannot be directly used as the clinical treatment basis.
- The electromagnetic field will affect the performance of this instrument, so the use of the other equipment near this instrument must meet corresponding EMC requirements. For example: Mobile phone, X-ray or MRI equipment may be an interference source, because they will transmit high-strength electromagnetic radiation.
- This is not a treatment device.

Caution

- To avoid damage to this instrument and guarantee patient safety, please use the accessories designated in this instruction manual.
- Please properly install or move this instrument and prevent the instrument from being damaged due to fall, collision, strong vibration or other external mechanical forces.
- Before the instrument is switched on, please confirm whether the power supply used meets the requirements for power supply voltage and frequency designated in the nameplate label or in the instruction manual of this instrument.
- When this instrument and its accessories are about to exceed the service life, they must be disposed of according to local relevant laws and regulations or the rules and regulations of the hospital.

 Attention

- Please install the equipment in a place that is convenient for observation, operation and maintenance.
- This instruction manual introduces the product according to the most complete configurations. The product you have purchased may not possess some configurations or functions.
- Please place this instruction manual near the instrument for easy and timely reference.
- This instrument cannot be used at home.
- This instrument can be used for one patient only at the same time.

Overview of the Monitor

This monitor (shown as follows) can be used for the bedside monitoring of adults, children and infants.

This monitor can monitor electrocardiogram (ECG), respiration (RESP), blood oxygen saturation (SpO₂), noninvasive blood pressure (NIBP), two-channel body temperature (TEMP), carbon dioxide (CO₂), anaesthetic gas (AG), noninvasive cardiac output (ICG), invasive blood pressure (IBP) or one of the main parameters.

The screen of this instrument adopts the 17-inch TFT color touch screen, and simultaneously supports two operation modes, i.e. key mode and touch mode. Next, we will introduce the basic functions of the monitor, shown as Figure 1-1:

④


Figure 1-1 Monitor

①	Display screen
②	Keys (from left to right: mute / reset key, freezing / unfreezing key, record / stop key, start/stop key for NIBP measurement, popup / exit the main menu key, flying shuttle mouse key). See the content of “Key Functions and Basic Operations” in this Chapter for more details.
③	From left to right: power switch, AC power indicator lamp, battery working status indicator lamp, battery charge indicator lamp.
④	Alarm indicator lamp

Screen Display

The screen of this monitor adopts the touch color LCD screen and simultaneously displays the patient parameters, waveforms, alarm information, clock, network connection status, bed number, battery and other prompt messages, etc.

The main screen is divided into four areas: 1. information prompt area or upper menu bar area ①; 2. parameter area ②; 3. lower menu bar area ③; 4. waveform area ④, which are shown as follows:


Figure 1-2 Standard Interface

Information Prompt Area(①):

The information prompt area is divided as follows from left to right:

- (a) Network connection state and network bed number: connection state of central monitoring system: when the symbol “” is displayed, it indicates that the central monitoring system is not connected; when the symbol “” is displayed, it indicates that the central monitoring system has been successfully connected and the figure next to it is network bed number. Click here to rapidly enter into the menu “Network Bed Number”.
 - (b) Information prompt of patient type (adult, child or infant). Click here to rapidly enter into the menu of “patient information” and set the detailed information about the patient. See the content of “3.1 Patient Information” for more details.
 - (c) When the pacemaker is turned on, it will be indicated by “” above the electrocardiogram waveform, and the symbol “” will be displayed at the top right corner of the information prompt area; when the pacemaker is turned off, there is no information prompt.
 - (d) Clock: display the current system time of the monitor. Click here to rapidly enter into the menu “System Time Setting” and then reset the system time of the monitor according to the local time zone.
 - (e) Alarm volume off: when the “alarm volume” is turned off and set as Level 0, the symbol “” will be displayed; when the “alarm volume” is set as other levels, there is no information prompt.
- Mute prompt. If the mute key is pressed down now, the prompt is usually: permanent mute or alarm suspending for 119s. The prompted content shall be determined by the setting of “Alarm Suspending Time” in the “Alarm Setting”, with the symbol “” displayed aside. If there is the symbol “” now, the symbol “” will be shielded.
- (f) Check the technological alarm, e.g. shedding of ECG lead and no connection to the probe for SpO2. Click here to rapidly enter into the window “Alarm Event Review” of technological alarm and then set the starting time of alarm in the window, and check the alarm information by pulling up or down the menu.

(g) Information prompt or mute prompt of physiological alarm

Information prompt of physiological alarm, e.g. ***RR too high and ***HR too low. Click here to rapidly enter into the window “Alarm Event Review” of physiological alarm and then set the starting time of alarm and the alarm event to be displayed in the window. For example, the events which can be selected include: All, ECG, SpO₂, NIBP, CO₂, RESP or TEMP, etc. Check the alarm information by pulling up or down the menu.

Parameter Area (②):

The parameter displayed in the parameter area is not fixed. The parameter on the right is displayed according to the waveform of the left. Usually, the waveform displayed on the left channel and the parameter on the right are displayed on the same channel, and the waveform and the parameter color are the same.

Lower Menu Bar Area (③):

Lower menu bar area is divided from top to bottom and from left to right: information prompt area, SD card status area, battery icon, mute key, alarm pause key, screen highlight menu, alarm volume menu, heartbeat volume menu, patient management menu, interface switch menu, review menu, alarm setting menu, NIBP measurement key, freezing key, list interface, respiratory oxygenation chart interface tendency coexistence interface, full screen 7-lead interface, full screen 12-lead interface, print key, default setting, monitor maintenance, measurement setting and main menu, etc. A part of them belongs to shortcut keys, which will not be described in detail here. Please refer to the relevant contents in each chapter.


Figure 1-3 Lower Menu Bar

Click the main menu, and the screen can be locked. When here is clicked and held for more than three seconds and then the screen changes from “” to “”, it indicates that the screen has been locked. In case of unlocking, continuously click the symbol “” for more than three seconds. When the symbol turns from “” to “”, the locked screen can be unlocked.

Waveform Area (④):

The waveform displayed in the waveform area is not fixed, which is related to the module equipped with the instrument. The more the monitoring parameters equipped with the module, the more the waveforms displayed in the waveform area. In the normal interface, if all modules are equipped with the instrument, at most the waveforms of 8 channels can be displayed.

The name of waveform on each channel is displayed on the top left of the waveform. ECG waveform can select the gain and filter modes to be displayed as required. For detailed information, please refer to the chapter “ECG Monitoring”. Except ECG waveform, other waveforms do not display the currently selected gain. There is a 1-mV scale rod at the right side of ECG waveform. Diminished or magnified multiple of scale rod varies with the gain.

During the screen operation, the menu can be popped up at a fixed location and the menu will occupy a small part of the location of the whole interface and thus part of the waveforms can be temporarily invisible.

The display of the original picture can be recovered after menu exit.

The waveforms shall be refreshed at a set rate. See the setting content of each parameter for the adjustment of refreshing rate of each waveform.

Key Functions and Basic Operations

Key Functions

The operations on this monitor can be finished by keys and knobs, shown as the following table:

Key symbol	Description
	(Mute/reset key) Alarm suspending for 10 minutes (this option can be set in the menu “Alarm Setting” in the main menu). Press down this key for more than one second to shield or turn off all alarm sounds, including the technological alarm, and display “Alarm Suspending for XXs” or “  ” in the information area. Under the status of “Alarm Suspending” or “Permanent Mute”, if a new technological alarm occurs, the current status of alarm suspending or permanent mute will end and the alarm will be recovered. After the alarm resumes, the instruct will immediately respond to the alarm and the response time will not exceed 3s.
	(Freezing/unfreezing key) In normal mode, press down this key to freeze all the waveforms on the screen.Press down this key again to release the frozen waveforms. After pressing down the freezing key, click the screen to review the waveforms stored within 287s.
	(Record/stop key) Press this key to start one real-time recording. The record length is determined by the “Real-time Recording Time” of the menu “Record Output Setting” in the “Main Menu”. The real-time recording time is respectively: 3s, 5s, 8s and continuous. During the continuous recording of recorder, press this key again to stop recording immediately. The key pressing frequency shall not be too high, and the interval of key pressing shall be greater than 2s.
	(Start/Stop Key for NIBP Measurement) At the non-blood pressure measurement state, press this key to aerate the cuff and then start one blood pressure measurement; at the measurement state, if you want to quit such measurement, press this key to stop measurement and then deflate.
	(Pop up/exit the main menu) In any interface mode, press this key to pop up or exit the main menu.
	(ON/OFF Key) Press this key to execute on-off control for the instrument.

Basic Operations

(1) Display the required waveforms:

In the normal interface, if all modules are equipped with the instrument, at most the waveforms of 8

channels can be displayed. However, if the module you purchase is only a certain one, the waveform and parameter displayed in the monitor interface are only the functions of a certain part.

Click the waveform on the channel you want to exchange. In the popup menu, select the menu “Exchange Waveform” to select the waveform you want to exchange.

(2) Adjust the waveform speed:

Except the two-channel ECG waveform, the waveform speed on each channel can be independently set, so you can select the appropriate scanning speed as needed. Click the waveform on the channel you want to exchange. In the popup menu, select the menu “Waveform Speed” to select the waveform speed you want to exchange.

(3) Change the alarm limit:

The alarm limit setting of each parameter can be set via its menu. Click the content to be set in the screen parameter area of the waveform area, enter into the menu, and select the menu of alarm limit. If the symbol “” is displayed on the parameter interface, it indicates that the alarm setting is located at the “Unconnected” position. Click the “Connected / Unconnected” switch to set as “Connected”, and set the alarm values of upper and lower limits.

(4) Adjust the volume:

The alarm volume, heartbeat volume and key volume can be set via the corresponding menus of the lower menu bar of interface. For example, click the “Alarm Volume” of the lower menu bar, and set the volume as Level “1~6”. The greater the value is, the greater the sound is. Vice versa. “0” indicates that the volume is “OFF”. Other settings may be deduced by analogy.

(5) Set the system time:

Click the clock area in the upper menu bar of the screen, and then set the time in the menu “System Time Setting” popped up.

External Interface of the Monitor

Left Panel of the Monitor

The following interfaces are provided on the left panel of the monitor, shown as follows:


Figure 1-4 Left Panel

! Attention

- Only the C60 products produced by the Company can be inserted into the small host slot.
- This monitor only can be compatible with the modules produced by the Company.

Right Panel of the Monitor

The following interfaces are provided on the right panel of the monitor, shown as follows:


Figure 1-5 Right Panel

Warning

- This network port only can be connected with the central monitoring system of Comen Company.
- All the simulation and digital equipment connected with this monitor must be the products certified by the designated IEC standards (e.g. IEC 60950 Data Processing Equipment Standard and IEC 60601-1 Medical Equipment Standard). Moreover, all configurations shall abide by the content of the valid edition of IEC 60601-1-1 System Standard. Connect the additional equipment to the staffing medical system at the input/output signal port and confirm whether the system conforms to the IEC 60601-1-1 Standard. If you have any question, please contact the supplier.
- When the signal interfaces like patient cable interface and network interface are simultaneously connected with multiple equipments, the total leakage caused cannot exceed the tolerance.

Back Cover of the Instrument


Figure 1-6 Back Cover

Equipment Symbols

(1) Instrument Symbols

	Attention! please see the accompanying documents	 xx-xx	Type approval mark and number of measuring instrument
	The application part of Type CF has the anti-defibrillation function	 Y.Z.No.:xx	Production license mark and number of measuring instrument
	The application part of Type CF		AC indicator lamp
	The application part of Type BF		Production Date Mark
	Start/Stop Key		Serial number mark
	Battery working state indicator lamp		Equipotential symbol
	Battery charging indicator lamp		Network connection symbol
	DVI interface		Multifunctional socket
ECG Defib 	Synchronous interface for defibrillation	Nurse Call	Nurse calling socket
	USB interface		SD card interface
	Offset calibration (reserved)		Main menu (reserved)

Note: Please see the content of “1.4 Key Function and Basic Operation” for the key symbols and their functions of the monitor.

(2) Packaging Symbols

	Up		Limit of stacking layers
	Fragile		Rainproof

Chapter 2 Installation of the Monitor

Attention

- To ensure that the monitor works properly, please read the information in this chapter, security information and patient safety chapters before using, and install monitor as required.

2.1 Unpacking and Examination

Carefully unload monitor and accessories from the box, and save the packaging materials for later transport or storage. Please compare the accessories according to the packing list. Check to see if there is any mechanical damage. Check all the external wires, insert any accessories needed. If there are any questions, please contact our sales department or agency immediately.

2.2 Connect the AC power cord Steps of connection

Confirm that the AC power supply meets the following specifications: a.c. 100V~250V, 50Hz/60Hz.

Use the power cord and power adapter equipped with the monitor. Insert the power cord into the interface of power adapter. Connect the other end of power adapter to the interface of monitor DC15V, and then insert the other end of power cord into the grounded three-phase power socket.

Attention

- Plug the power cord to the dedicated hospital outlet.
- If a battery is provided, you must charge the battery after the transport or storage of instrument. If you do not connect the AC power and directly turn on the monitor, it will probably not work because of insufficient battery power. Connect to an AC power supply and you can charge the battery ,regardless of whether the monitor is turned on or not.

Connect the equipotential grounding wire when necessary. Refer to the content of equipotential grounding in the chapter “Patient Safety”.

2.3 Power on After you turn on the power switch, the screen will display “system is loading ...”.

then it will display the Company's LOGO, perhaps for 1-5 seconds. After the systems self-test has succeeded and entered the main screen, user can operate at this time.

 Attention

- If there are significant errors in the self-test process, the system will alert.
- Check all monitoring functions that can be used to ensure that the monitor is functional.
- If a built-in battery is inside the unit, the battery must be recharged after each use to ensure enough energy in the battery.
- Power on again after 1 minute of power off.

 Warning

- If you find signs of damage to monitor functions, or an error message, do not use this monitor for patient monitoring. Please contact the biomedical engineers of your hospital or maintenance engineer of the Company.

2.4 Connecting the Sensor Connect the required sensor to the monitor and the monitoring part of patient.

 Attention

For the correct connection methods for sensors and related requests, see Chapters 10-15.

Chapter 3 System Settings

3.1 Measuring Settings

You can select the “Measuring Settings” in the lower menu bar area through the touch mode or through pressing the key “↩”, shown as in the figure below:


Figure 3-1 Measuring Settings

What is displayed on this menu is the module information of the current plug-in unit. For example, the modules used currently include AG module, CO2 module and C30-type monitor.

- Disabling of certain module or certain monitoring function: When certain module or certain monitoring function is not used, set as “Disabled” status, with the steps as follows:
 - (a) Enter Shortcut, remove “√” for the module or certain monitoring function to be disabled, and select Yes in the popped-up dialog box Prompt.
 - (b) Click the module or certain monitoring function to be disabled, and select Yes in the popped-up dialog box Prompt.
 - (c) “×” will be displayed on the disabled module or certain monitoring function.
- If you want to recover the use of certain module or certain monitoring function, follow the above reverse steps.

3.2 Main Menu The system settings for this monitor are very flexible. You can enter the window “Main Menu” through the touch mode or through pressing the key “↶”, and then set each menu of the system via this window, shown as in the figure below:


Figure 3-2 Main Menu

 Attention

- The system settings of the monitor, once set, will be storage and memorized, and will not change because of power failure or instrument closing unless artificial or re-default setting.

3.2.1 Patient Management

3.2.1.1 Receiving (Admitting) Patients

 Attention

- When “admission” of a new patient, please confirm whether the original patient has been discharged or changed to another department, and correctly remove (discharge) such patient; otherwise, the data of the next patient will be stored in the data of such patient.
- After turning on the “Pacemaker”, the system will always keep the last settings unless “Factory Default Setting” is carried out again.

When connecting a patient with the monitor, even if no patient is received (admitted), the monitor can also

display the stored physiological data of the patient, but it is of great importance to correctly receive (admit) patients.

Receiving an admitted patient:

- (a) Press the key “↩” to enter Main Menu, and select Patient Management → Admission.
- (b) Select YES or NO to remove the current patient and receive a new patient.
 - ✧ YES: A new patient has been successfully received, and all the measuring results and historical data of the original monitored patient have been removed.
 - ✧ NO: Exit the receiving of a new patient.
- (c) Input the information of patient in the menu Patient Information, and in particularly, correctly select Patient Type and Pacemaker:
 - ✧ Patient Type: It is rather important to correctly select the option, which determines the measuring algorithm of calculation and processing used by the monitor as well as the scope of certain safety limits and alarm limits.
 - ✧ Pacemaker: When the pacemaker is turned on, it will be indicated by “I” above the electrocardiogram waveform, and the symbol “❤️” will be displayed at the top right corner of the information prompt area; when the pacemaker is turned off, there is no information or symbol prompt.
- (d) The menu is as shown in the figure below:

The screenshot shows a menu titled "PAT MANAGE" with a close button (X) in the top right corner. The menu contains the following fields and values:

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At the bottom of the menu, there are two buttons with left and right arrow icons.

Figure 3-3 Patient Information

 Warning

- When the patient type changes, the upper and lower alarm limits of some parameters will change accordingly. Usually, the alarm limits should be verified before patient monitoring to ensure that these alarm limits suit your patient.
- For non-pacing patients, you should verify that setting is "off"; otherwise, the system cannot detect arrhythmias associated with ventricular premature beats (including PVCs count) as well as ST segment analysis.
- For pacing patients, pacing pulse analysis features must be started. Otherwise, pacing pulse may be counted into normal QRS waves, causing an alarm of "ECG signal is too weak" and cannot be detected.

 Attention

- When patients are monitored with pacemakers, "pacemaker" (PACE) switch must be turned on; When patients are monitored without pacemakers, "pacemaker" (PACE) switch should be off; When "pacemaker" switch is on, some kind of ARR analysis will not be conducted on by the monitor. For more information please read the introduction section of the arrhythmia chapter.
- When pacing analysis is on, arrhythmias associated with premature ventricular beats (including PVCs count) will not be detected, as with ST segment analysis.

3.2.1.2 Removing (Discharging) Patients

Removing a discharged patient:

- (a) Press the key  to enter Main Menu, and select Patient Management → Discharge.
- (b) Select YES or NO to remove the current patient.
 - ✧ YES: The current patient has been successfully removed, and all the measuring results and historical data of the original monitored patient have been removed.
 - ✧ NO: Exit the removing of the current patient.

3.2.2 Alarm Settings

You can enter the window “Main Menu” through the touch mode or through pressing the key “↩”, and then select the menu Alarm Settings in this window for alarm settings. For more information on alarm settings, please refer to the chapter “Alarm”.

3.2.3 Recording Output Settings

Recording output settings have the following options:

- Record the waveform on the channel 1 to select any one ECG lead you want to print or the waveforms of other monitoring functions, e.g. Pleth, Resp, CO₂, O₂, N₂O and AA, etc..

- Record the waveform on the channel 2 to select any one ECG lead you want to print.

- Recording output speed:

There are two options: 25 m/s and 50mm/s.

- Recording network:

It is used to set the record output format. If the setting “OFF” is selected, it indicates that there is no network output; if the setting “ON” is selected, it indicates that there is network output.

- Real-time recording time:

There are four options: 3s, 5s, 8s and Continuous. Continuous: Once the button “” on the panel is pressed, the recorder will continuously output waveforms and parameter values until this button is pressed again.

Attention

- The recorder is an optional component.

3.2.4 Review

You can enter the window Main Menu and select the menu Review through the touch mode or through pressing the key “↩” to review the tendency table, tendency chart, NIBP, alarm event and waveform holography. Please refer to the chapter “Review” for details.

3.2.5 Drug Calculation

The drug calculation and titration meter function of the monitor can provide the drug concentration calculation of 15 different drugs. For the application method, refer to the chapter “Drug Calculation and Titration Meter”.

3.2.6 Interface Switch

The monitor provides several interfaces for the user, e.g. normal (standard) interface, 4-waveform interface, 6-waveform interface, 8-waveform interface, large font interface, tendency coexistence interface, respiratory oxygenation chart interface, list interface, full screen 7-lead interface, full screen 12-lead interface and other bed observation interface. One option can only be selected once. Please refer to the chapter “System Working Interface” for details.

3.2.7 Monitor Information

Monitor information is mainly for checking the version information of the software operated by the instrument, so as to facilitate the manufacturer to maintain and trace the instrument.

3.2.8 Monitor Maintenance (user maintenance)

You can enter the window “Main Menu” through the touch mode or through pressing the key “”, and then select the menu Monitor Maintenance. In the menu Monitor Maintenance, enter the correct password: 5188, and press the button “Enter” to enter the maintenance menu interface.

- Language Selection: Chinese/ English.
- Monitor Resolution: Default 17 TFT.
- Touch screen calibration: select "touch calibration" menu. The screen calibration interface appears. According to the screen calibration point that pops up “”, as the confirmation that five calibration points “” have been clicked. Calibration ends and the screen automatically returns to the main interface. Nurse calling setting:

Alarm Switch: “ON” / “OFF”.

Alarm Level: “High”, “Medium” and “Low”.

- Factory default setting: select this item and enter "factory defaults" dialog box. Selecting "no" means the current operation is aborted and the original configuration of the system is restored. Selecting "Yes"

means the factory default settings of adult configuration is to be used and the original configuration will be overwritten.

- After the factory setting is overridden by default setting, all the monitoring parameter alarm preset values will be redefined. Therefore, when using this function, please check whether the type of patient you are monitoring currently is correct.

Attention

- Default setting means that if you choose "Default Setting", the system will override the current parameter with the default parameters provided by manufacturers.
- The user cannot execute the factory maintenance setting function in the monitor maintenance. This option is opened only for the maintenance personnel designated by the Company.

3.2.9 Demo Features

Select "Demo features" in the menu "Monitor Settings" and a pop-up window "Enter Your Password" dialog box will appear. After entering the correct password, select "ON", and then the monitor enters presentation status.

Warning

- Demo waveform is a kind of simulation of waveforms, which is made by the manufacturer only to demonstrate the machine performance, and help users to set up training. In actual clinical use, the presentation is disabled, because it may make the medical staff mistake it for the monitoring waveform and parameters of patient monitored, which affects monitoring and delays diagnosis and treatment. So this menu has a password.

3.2.10 System Time Settings

Push " key to enter "shortcut key" window, select "main menu" in this window and enter "system time setting" menu or click the clock area in the menu bar below the screen in the interface. In the pop-up "system time settings" menu, set system time in accordance with the local time zone. The time settings include: year, month, day, hour, minute and second.

3.2.11 Network Bed Number

Push “↩” key to enter the "shortcut key" window, select "main menu" in this window and enter "network" menu. From this menu, you can set up network bed number, IP address, MAC address, subnet mask, server IP, etc, which is connected to the central monitoring system of the company. In General, we only need to set the network number, and other settings are in the default state.

- Network bed number: bed number of monitors connected to the central monitoring system network.
- IP address: 200.200.200.X, (x, refers to the network bed number, 1~128 beds are optional).

How to know the network is connected successfully? There is a central monitoring system icon below the interface area. When it appears “”, it indicates that the central monitoring system has not connected successfully. When it displays “”, it indicates that the central system has connected successfully, next to the icon is the network number.

Attention

- The Network bed number must be unique, and can not overlap the bed numbers of any other instruments connected to the central monitoring system, otherwise it will cause instrument signal deadlock because of the preemption of the central monitoring system channel.
- If any instrument freezes due to network bed number conflicts, remove the network cable, turn off the monitor and restart. Reset the networks and then reconnect to the network connector.

Chapter 4 System Working Interface

4.1 System Working Interface

The monitor provides several interfaces for the user.e.g. normal(standard) interface, 4-waveform interface, 6-waveform interface, 8-waveform interface, large font interface, tendency coexistence interface, respiratory oxygenation chart interface, list interface, full screen 7-lead interface, full screen 12-lead interface and other bed observation interface. The user can get different screen messages by choosing a different work interface based on different requirements. Next is the introduction to the styles and characteristics of part of working interfaces.

4.1.1 Normal Interface

Normal interface is also known as standard interface. If the instrument has complete functions, at most waveforms of 8 channels can be displayed on this interface, with divided parameter display area. The waveform on each channel is not fixed. The waveform displayed on the interface is displayed according to the equipped function of your instrument, and you can select the channel where the waveform displays as needed.

Open the normal interface:

- (a) In the lower menu bar area, or enter the Main Menu to select Interface Switch;
- (b) To select the Normal Interface, as shown in the figure below:


Figure 4-1 Standard Interface

4.1.2 Large Font Interface

Large font interface can display four parameters and four waveforms. The measured value is displayed with super-large font to facilitate the doctor to observe the monitored value within a certain distance.

Open the large font interface:

- (a) In the lower menu bar area, or enter the Main Menu to select Interface Switch;
- (b) To select the Large Font Interface, as shown in the figure below:


Figure 4-2 Large Font Interface

4.1.3 Tendency Coexistence Interface

What is displayed on the tendency coexistence interface is the dynamic short tendency chart of waveform. From the dynamic short tendency chart, you can see the recent tendency of a series of parameters. The mark of the tendency is displayed above the tendency chart, with the scale of the tendency displayed on the left, and the time of the tendency chart displayed below.

The dynamic short tendency chart position of each parameter is not fixed. You can select the appropriate tendency chart as needed. Take the HR dynamic short tendency for example, with the steps as follows:

Select the HR dynamic short tendency (see the figure below), enter the menu Tendency Parameter Setting, and select the SPO₂. Exchange for SPO₂ dynamic short tendency at the position of the original HR dynamic short tendency.


Figure 4-3 HR Dynamic Short Tendency

Open the dynamic short tendency interface:

- (a) In the lower menu bar area, or enter the Main Menu to select Interface Switch;
- (b) To select the Tendency Coexistence, as shown in the figure below:


Figure 4-4 Tendency Coexistence Interface

■ Position of tendency chart

Tendency chart is located on the left of the corresponding waveform in the waveform area, whose color is the same as the corresponding parameter color.

■ Tendency length

The dynamic tendency length is 2h, with the right end of abscissa of tendency chart being 0h, and the left end being 2h.

■ End the tendency coexistence interface

In the selection menu Interface Switch, select other working interfaces to end the tendency coexistence interface.

4.1.4 Respiratory Oxygenation Chart Interface

Respiratory oxygenation chart occupies one half of the entire waveform area, which is comprised of HR tendency, SPO2 tendency and RR tendency or compressional respiratory wave.

Open the respiratory oxygenation chart interface:

- In the lower menu bar area, or enter the Main Menu to select Interface Switch;
- To select the Respiratory Oxygenation Chart, as shown in the figure below:


Figure 4-5 Respiratory Oxygenation Chart Interface

■ Tendency chart of respiratory oxygenation chart interface

Respiratory oxygenation chart interface is comprised of two tendencies, i.e. “compressional respiratory wave” and “RR tendency”, located below the screen waveform area.

■ Respiratory oxygenation chart tendency time span selection

Select the respiratory oxygenation chart, and one menu Respiratory Oxygenation Chart will pop up:

- ✧ Time: You can select the time spans: 1min, 2min or 4min.

■ Compressional respiratory wave and RR tendency selection

Select the respiratory oxygenation chart, and one menu Respiratory Oxygenation Chart will pop up:

- ✧ Type: Compressional respiratory wave or RR tendency

The tendency function of “compressional respiratory wave” is according to the different display modes of “RR tendency” or “compressional respiratory wave” selected by the operator. The displayed contents under the two modes occupy the same position. When the “RR tendency” is selected, the position will display the dynamic tendency chart of respiratory rate; when the “compressional respiratory wave” is selected, the position will display the compressed respiratory wave.

- End the respiratory oxygenation chart interface:

In the selection menu Interface Switch, select other working interfaces to end the respiratory oxygenation chart interface.

4.1.5 List Interface

Measurement parameter list occupies one half of the entire waveform area. On this interface, you can review the recently measured data, which is stored based on 12-group data per page.

Open the list interface:

- (a) In the lower menu bar area, or enter the Main Menu to select Interface Switch;
- (b) To select the List Interface, as shown in the figure below:


Figure 4-6 List Interface

4.1.6 Other Bed Observation

Other bed observation means that the monitor can observe the information of any measurement waveform and all measurement parameters of another monitor connected to the same network. Other bed observation

window will occupy 4-channel waveform display area at the bottom of screen.

 Attention

- **This function is optional, and can only observe the machine of Comen Company.**

Open the other bed observation interface:

- (a) In the lower menu bar area, or enter the Main Menu to select Interface Switch;
- (b) To select the Other Bed Observation, as shown in the figure below:


Figure 4-7 Other Bed Observation Interface

4.1.7 Other Interfaces

4-waveform interface, 6-waveform interface and 8-waveform interface can display the waveforms of 4 channels, 6 channels and 8 channels respectively. The less the number of channels is, the larger the distributed area of waveform is, and you can observe more clearly.

Full screen 7-lead interface can display the ECG waveforms of 7 leads “I, II, III, AVR, AVL, AVF, V”, applicable to the comprehensive observation of ECG waveform when the lead type is “5-lead”.

Full screen 12-lead interface can display the ECG waveforms of 12 leads “I, II, III, AVR, AVL, AVF, V1, V2, V3, V4, V5, V6”, applicable to the comprehensive observation of ECG waveform when the lead type is “12-lead”.

4.2 Setting the Interface Style

You can set the interface style as needed, including:

- ◆ Waveform scanning mode;
- ◆ Waveform display pattern;
- ◆ Parameter/waveform display color;
- ◆ Waveform to be exchanged;
- ◆ Gain.

The setting method of each parameter of the monitor interface style is the same. In this chapter, only the setting method of one parameter will be illustrated, and the setting methods of other parameters will not be introduced repeatedly.

4.2.1 Setting the waveform scanning speed

- (a) On the normal interface, click the waveform you want to set. For example, select the waveform II or I of ECG;
- (b) In the popped-out menu ECG Waveform_i, select the Waveform Speed;
- (c) Set the Waveform Speed as 6.25mm/s, 12.5mm/s, 25.0mm/s or 50.0mm/s;
- (d) Upon completion of setting, press the key “X” to exit.

4.2.2 Setting the waveform display pattern

- (a) On the normal interface, click the waveform you want to set. For example, select the waveform II or I of ECG;
- (b) In the popped-out menu ECG Waveform, select the Waveform Pattern;
- (c) Set the Waveform Pattern as Line, Level or Bold;
- (d) Upon completion of setting, press the key “X” to exit.

4.2.3 Setting the display color

- (a) On the normal interface, click the waveform you want to set. For example, select the waveform II or I of ECG;
- (b) In the popped-out menu ECG Waveform, select the Waveform Color;
- (c) Set the Waveform Color as Green, Cyan, Red, Yellow, White, Blue, Violet or Orange;
- (d) Upon completion of setting, press the key “” to exit.

4.2.4 Setting the waveform to be exchanged

- (a) On the normal interface, click the waveform you want to set. For example, select the waveform II or I of ECG;
- (b) In the popped-out menu ECG Waveform, select the Exchange Waveform;
- (c) Set the Exchange Waveform as the waveform you want to exchange to display parameters;
- (d) Upon completion of setting, press the key “” to exit.

5.1 Alarm Type

The monitor can give an alarm of three types: physiological alarm, technical alarm and prompt message.

(1) Physiological alarm

Physiological alarm is usually caused by a certain physiological parameter of the patient which exceeds the set upper/lower limit scope or by the physiological abnormality of the patient. The alarm message of the physiological alarm will appear in the physiological alarm area above the screen.

(2) Technical alarm

The technical alarm is also referred to as a system error message, indicating the alarm is caused by a misoperation or system malfunction thereby causing improper operation of a system function or distortion of monitored results. The alarm message of the technical alarm will appear in the technical alarm area above the screen.

(3) Prompt message

Strictly speaking, prompt message does not belong to the domain of alarm. It is to display the information relative to the system conditions themselves, which have nothing to do with the patient's vital sign. Prompt message will appear in the prompt message area.

5.2 Alarm Level

The monitor can give an alarm in three volume levels: high, medium and low, based on the severity.

(1) High level alarm:

The patient is in a critical condition, endangering the patient's life, emergent attention required.

(2) Medium level alarm:

The patient's vital signs are abnormal, relevant measures and treatment are immediately required.

(3) Low level alarm :

The patient's vital signs are abnormal, and relevant measures and treatment may be required.

All alarm levels for technical alarms and some physiological alarms have be set before the monitors

are delivered:, the users are not allowed to change them. but a certain degree of physiological alarm can be modified.

5.3 Alarm Mode

When giving an alarm, the monitor will prompt the users in both sound and visual modes:

Lighting Alarm

Sound Alarm

Alarm

parameter's flashing

In which light signal, sound signal and prompting message are differentiated with different alarm levels.

5.4 Lighting Alarm When giving an alarm, the alarm indicating lamp will show alarms in different levels with color and flashing frequency.

High level alarm red with flash frequency of 2 per second.

Medium level alarm yellow with flash frequency of one every two seconds.


Low level alarm yellow without flash but with a continuous light.

5.5 Sound Alarm The sound alarm is set to prompt the alarm in different levels of severity with different sounds.

High level alarm: Do-Do-Do--Do-Do----Do-Do-Do--Do-Do, response once at the interval of 6s.

Medium level alarm: Do-Do-Do, flash once at the interval of 10s.

Low level alarm: Do, flash once at the interval of 10s.

- Both bedside machines and central monitoring systems have sound alarming function.
- Once the bedside machine is connected to a central monitoring system ,although both the bedside machine and central monitoring system can coordinate the alarm level as well as the upper and lower limits of the alarm, the bedside machine may not give an alarm simultaneously when the central monitoring system is prompting an alarm due to the delaying function of the bedside machine.

5.6 Alarm

Prompting message means the related information will appear on the physiological alarm area or technical alarm area of the monitor when the alarm is active. The system will show different alarm levels in different background colors:

High level alarm: Red

Medium level alarm: Yellow

Low level alarm: Yellow

Levels of the prompt message in the front of the alarm will be distinguished with the following symbols:

High level alarm: ***

Medium level alarm: **

Low level alarm: *

5.7 Alarm Parameter Flashing When a parameter is alarming, the parameter will flash once every second.

5.8 Icon of Alarm Status Besides the alarm modes mentioned above, the following alarm icons will appear on the screen to indicate different alarm conditions.

: indicating alarm muted.

: indicating alarm pause.

: indicating shut off of a parameter alarm.

: indicating Level 0 OFF mode of alarm volume.

5.9 Setting of Alarm Press “↶” key to go to the “shortcut key” window to choose the main menu, then go to the menu “Alarm Setting” as shown below:


Figure 5-1 Alarm Setting

- Alarming volume: five alarming volumes are available in the system: “turn off”, and 1-4 levels., Once the 1-4 levels selected the system will prompt the alarm at a certain alarm volume, Accordingly; while “turn off” is selected will shut off the alarm sound.
- In “time-out of alarm”, there are: six different timeouts provided with the system: “permanent”, “one minute”. “two minutes”, “three minutes”, “five minutes” and “ten minutes”. This function is effective only when the system is giving an alarm, but when a new technical alarm is given the instrument will immediately prompt with an alarm, as if the “time-out of alarm” is disfunctional. When the “permanent” is selected, when the “sound off” key is pressed, it prompts the permanent sound off. When the “one minute” time-out key is selected, press the sound off key to prompt a 59 second time-out for alarm, and make a countdown. and so on.
- Alarm delay: For the alarm delay time of the parameter, the system provides five options: Disable, 5 Seconds, 10 Seconds, 15 Seconds, and 20 Seconds. When “Disable” is selected, the instrument will give the alarm when the measurement value exceeds the alarm limit. When “5 Seconds” is selected, the instrument will give the alarm only when the measurement value exceeds the alarm limit for 5 seconds.

5.10 Set the Alarm Volume(a) Select the Alarm Volume in the lower menu bar area, select the volume of Level 0~6, and select OK.

When the alarm volume is set as 0, an icon  appears on the information prompt area of the screen, indicating that the volume is OFF.

 Warning

- When the alarm volume of system is set as “0”, the monitor will not make any alarm sound even if a new alarm event happens. Therefore, the operator should use the function carefully.

5.11 Setting parameter alarm

5.11.1 Setting the upper and lower limits of alarm as well as the ON and OFF alarms

Take ECG for example:

- Select ECG in the parameter area;
- In the popped-out menu ECG Setting, select the Alarm Limit;
- Set the upper and lower limits of alarm as 100bpm and 60bpm respectively;
- Set the Connected / Unconnected as ON.
 - ON indicates that when this parameter gives an alarm, the monitor will prompt according to the set alarm level, and store the relevant waveform and value of this parameter. At the same time, the upper and lower limits of alarm will display in the parameter area.
 - OFF indicates that when the alarm scope of this parameter exceeds the limit, the instrument has no alarm prompt or storage because the alarm switch is not turned on. There is an icon  prompt at the setting of Alarm Limit of this parameter and in the parameter area. At the same time, the upper and lower limits of alarm are turned off.
- Upon completion of setting, press the key “” to exit.

 Attention

- It is of great importance to set the upper limit of heart rate alarm during the monitoring. The upper limit should not be set too high. With a view to the fluctuating factors, the set upper limit of heart rate alarm is not 20bpm higher than the patient's heart rate.

Warning

- When setting the upper and lower limits of alarm, decide whether the patient you are monitoring is an adult, child or newborn, and the limit should be set according to clinical requirements. If the set limit value exceeds the alarm limit, it will easily lead to alarm system failure.
- When the alarm limit is turned on, after the upper and lower limits of alarm are set manually, the instrument will display the set upper and lower limits of alarm continuously, and will not additionally provide the initial alarm pre-set value of the system.

5.11.2 Setting the alarm level

Take ECG for example:

- Select ECG in the parameter area;
- In the popped-out menu ECG Setting, select the Alarm Level;
- Set Alarm Level as High, Medium or Low;
- Upon completion of setting, press the key “” to exit.

5.11.3 Alarm Record

If the monitor is equipped with a recorder, when this parameter gives an alarm, as long as the Alarm Limit (Connected / Unconnected) and the Alarm Record are set as ON, the instrument will trigger the recorder to output the relevant waveform and value of this parameter. Take ECG for example:

- Select ECG in the parameter area;
- In the popped-out menu ECG Setting, select the Alarm Record;
- Set the Alarm Record as ON.
- Upon completion of setting, press the key “” to exit. Verifying the Alarm System Whether or not the

alarm system is working, can be detected through the light and sound alarm status indicator.

- (1) connecting a cable for blood oxygen measurement to the monitor;
- (2) go to “SpO₂ setting” to set the blood oxygen alarm in an “ON” position
- (3) setting “SpO₂ alarm at high and low limits” in 90% and 60%, respectively
- (4) go to “main menu” to select “setting alarm” in which the alarm volume” is set at one of the levels of “1-4”;
- (5) When the measured value exceeds the upper or lower limit of alarm, set the Alarm Level as High, Medium and Low to observe the changes in the sound, light and parameter flash of the instrument. See the contents of “Sound Alarm, Light Alarm, Alarm Information and Alarm Parameter Flash” in this chapter for details. Meanwhile, it indicates that SPO₂ is too high or too low in the physiological alarm information area.
- (6) If it warns under the normal circumstance which the parameter surpasses the maximum and minimum of alarm, the equipment will response with alarms of “light, voice, parameter flash and information” , and the average delay of alarm does not exceed 5 seconds.
- (7) Pull out the SPO₂ probe from the monitor, and the following information prompts out in the technological alarm information area: the SPO₂ probe comes off.

Attention

- When the alarm system gives several alarms at the same time, the system preferentially gives high level warning with light and voice.

5.13 Time-out or Disabled Alarm

When alarm is given, an operator can press “” key to suspend or shut off the alarming sound and suspending the alarm sound for all parameters being monitored and shut off technical alarms. The time-out will display on the upper-right corner of the screen, such as: “Alarm time-out: XXs”. The “Alarm setting” on main menu can be set, accordingly.

Warning

- When the system is set in shut off sound, the monitor will not make any alarm sound if an alarm event happens. Therefore, the operator should use the function carefully.

Chapter 6 Battery

6.1 General

A built-in rechargeable battery is installed in the medical monitor. When connected to the AC power supply, the battery will be automatically recharged, whether the unit is in an ON or OFF position, until the battery is fully charged. In case the power suddenly turns off, the system will be automatically powered on by the built-in battery thereby not interrupting the unit while working, and the indicator light for the battery will be lit after the power supply has been shut off for over 30 seconds.

The symbol “” will be shown at the bottom right corner indicating the power condition of the battery: the green indicates the power is still full or medium level, while yellow shows a low level; and red an extremely low level, alerting the user of the condition.

The battery icon shown on the screen indicates the current battery level:

Indicating fully charged battery.

Indicating the battery is charged, but not at a full level.

Indicating the battery is at a low level, requiring recharging the battery.

Indicating the battery is unavailable or malfunctioning..

Attention

- Please remove the battery if it will be out of service for a long time and store it properly.
- If a built-in battery is inside the unit, the battery must be recharged after each use to ensure enough energy in the battery.

Warning

- The battery liquid is harmful, In case the liquid contacts your skin or eye, wash it immediately with large amounts of clean water or seek medical advice.
- Keep the battery out of reach of children.

- The monitor is working it will automatically turn off power when power is too low, if the energy in the battery has nearly run out, the monitor will give a high level alarm with a continuous sound of “beep”. If the battery is at a low level and prompts a message at the message area of the screen such as “Low energy level of battery”. If this case happens, the battery should be immediately connected to AC supply to be recharged. If the battery is still being used, the monitor will automatically shut off power before the battery is fully used up (after alarming for one minute).

6.2 Installing Battery

Procedure of changing or installing the battery:

- (1) Turn off the monitor, and disconnect the power cord and other connection lines.
- (2) Place the monitor with the back up.
- (3) Unscrew the battery cover.
- (4) Take out the used battery and put the new one into the battery holder making sure the positive and negative poles match what is indicated on its shell.
- (5) Replace the holder and screw it on, and turn the monitor upright.

Warning

- Use only the supplier's designated battery.
- Don't dismantle the battery while the monitor is turned on.

6.3 Optimization and check of battery performance

Optimization of battery performance

When the battery is used for the first time at least two complete cycles of optimization of the battery should be carried out. A complete optimization cycle should be: uninterrupted charging battery until the power is full, followed by use until the battery is fully discharged and monitor is automatically shut off.

This will ensure the battery is in optimization process:

- (a) Disconnecting the monitor from the patient suspends all monitoring and measuring procedures.

- (b) The optimized battery should be kept in the battery holder of the unit.
- (c) When charging the battery, at least six hours of charging should be ensured until it is fully charged.
- (d) When you disconnect the AC power supply, the monitor is powered with the battery until the battery runs out and the monitor automatically shuts off.
- (e) This completes the battery optimization process.

(2) Check of Battery Performance

The service life of battery is changeable along with its storage, working environment charge cycles and service time. Even though battery is out of service its performance will gradually deteriorate.

Procedure for checking the battery is as follows:

- (a) Confirming whether or not the battery is damage. When the battery shows the symbol “”, it indicates the battery is damaged or not in the battery holder.
- (b) Checking whether the battery can be normally charged when the battery is connected to alternating current;
- (c) Disconnecting the monitor from the patient suspends all monitoring and measuring procedures.
- (d) When charging the battery, at least six hours of charging should be ensured until it is fully charged.
- (e) Disconnect the AC power supply, power on the monitor with the battery until it is fully discharged and the monitor shuts off automatically Record the start and stop time.
- (f) The period of battery discharge will reflect the battery performance.
- (g) Once the discharge period is down to 50% of the original time, it requires changing the battery.

- In order to extend the service life of the battery it is recommended to charge it every three months after a long dormant period so as to prevent overdischarge.
- Battery power supply loss depends on the configuration and operation of the monitor; for example, the unit will have a big loss of battery power if it is used to measure NIBP parameter often.

6.4 Battery Recovery If the battery shows apparent damage or is at an energy exhaustion condition, it should be exchanged immediately, and the old battery should be recovered and properly disposed of in accordance with relevant laws or rules and regulations for hospitals.

Warning

- Do not remove the battery or make it short-circuiting or put it into fire; otherwise, it would cause battery on fire, explosion, harmful gas leakage or other dangers.

7.2 W2

The monitor provides the 1200-hour tendency data, 5000 sets of NIBP measurement data, storage of 500 parameter alarm events, review of 100 times of arrhythmia, 40 sets of 12-lead ECG analysis results, and 24-hour waveform holography review for all monitoring parameters (the specific time is related to the stored waveform and its amount). Methods for observing these stored data are to be explained in this chapter.

Copyright © 2016 by John Wiley & Sons, Inc.

NIBP RETRO						X
	SYS	DIA	MAP	PR	TIME	
1	120	80	90	60	2012-05-14 16:02:28	
2	120	80	90	60	2012-05-14 09:26:05	

Figure 1-1 NIBP measurement review

At this window, you may review the data of the pulse rate.

7.2 Warning Event Review

This monitor can display the latest 500 parameter warning events in warning event review. Directly click the option Review in the lower menu bar, or press the key “↶” to enter the Main Menu, enter the menu Review, and then select Alarm Event Review, as shown in the figure below:


Figure 7-2 Warning Event Review

Users can define conditions for warning events review in this menu, which includes:

(1) Starting Time for Warning Review

Users can set the “starting time” of warning reviews in the option Starting Time; this set time is limited by the system if there is data beyond the set starting time.

(2) Event Selection of Warning Review

Users can select the parameters they want to review in the Event Type of Warning Review. Options available include “All” (warning events of all parameters), or any parameter selected from ECG, SpO₂, NIBP, CO₂, RESP, TEMP, IBP<1, 2>, AG, ICG and ARR.

(3) Warning Event Review

Displayed in the menu Warning Event Review is the information as follows:

Warning Occurrence Time (Format: Year-Month-Day: Hour: Minute: Second);

Event Types (Warning Level included);

7.3 Arrhythmia Review (no window display)

This monitor displays the 100 latest ARR reviews during arrhythmia review. Press “↩” to enter “Shortcut Key” window, select “Main Menu” in this window to enter the “Data Review” Menu, and then select “ARR Review”, specifically as shown in Figure below:


Figure 7-3 Arrhythmia Review

Users can set the conditions for arrhythmia review in this menu, which includes:

(1) Starting Time for Arrhythmia Review

Users can define the “starting time” of warning reviews in this option; this defined time is limited by the system if there is data beyond the defined starting time.

(2) Arrhythmia Review

Displayed in the menu Arrhythmia Review shows the information as follows:

Arrhythmia Occurrence Time (Format: Year-Month-Day Hour: Minute);

Event Types (Warning Level included);

Sequence No. (Format: x in y).

Arrhythmia events before or after are displayed on other pages and can be reviewed by pressing the button PageDown or PageUp.

 Attention

- In case that the arrhythmia events exceed 100, the monitor will only save the latest 100 events and the earlier ones will be discarded. A monitor with power-failure storage functions can store the latest 100 arrhythmia events in case of a power failure.

7.4Review of tendency chart

The tendency chart in the latest hour can be displayed with the resolution of one data every second or every 5 seconds; and the tendency chart in the latest 120 hours can be displayed with the resolution of one data in every minute, or every 5 or 10 minutes. Directly click the option Review in the lower menu bar, or press the key “↶” to enter the Main Menu, enter the menu Review, and then select Tendency Chart Review, as shown in the figure below:


Figure 7-4 Tendency Chart Review

The vertical axis represents measured value and the horizontal axis represents measured time. Described as follows are the displayed symbols used:

Symbol	Description
 	Page up or down to review the tendency diagram of other parameters not displayed

	on the current view.
 	Move the cursor one step left or right to review a tendency diagram by the time line of tendency database.
 	Move the cursor one page left or right to review a tendency diagram by the time line of tendency database.
 	Skip to the starting or ending point of a tendency database to review the furthest (earliest) or the nearest (latest) stored tendency information.

(a) Select and Display Tendency Diagrams of Different Parameters:

Press the key “” or “” to select the desirable parameter, and then the tendency chart of this parameter will appear in the window.

(b) Select 1h or 120h Tendency Chart:

Click on the option “Resolution”, and select 1 or 5 seconds to review a tendency chart for the last hour and 1, 5 or 10 minutes to review a tendency chart in the last 120 hours.

(c) Review Previous or Nearer Tendency Curves

Press “” to enter window on next page, and then press “” “”. Press the left button to review tendency curves on earlier time points, and the right button to review tendency curves on nearer time points.

(d) Obtain Tendency Data on a Given Moment from Current Diagram

Press “” to enter window on next page, and then press “” “” to move cursor, and the selected moment will change when the cursor moves, the values corresponding to the parameters at the selected moment will be displayed on the horizontal axis, and the values on the horizontal axis will be displayed on the right of the window.

(e) Operation Example

Review NIBP Tendency Diagram in the Last Hour:

- (1) Directly click the option Review in the lower menu bar, or press the key “” to enter the Main Menu, enter the menu Review, and then select Tendency Chart Review;
- (2) In the window Tendency Chart Review, press the key “” or “” to select the desirable parameter;

- (3) Select “1 sec.” or “5 sec.” in the “Resolution” Option.
- (4) Press “▶▶” to enter window on next page, and then press “◀◀” “▶▶” or “◀” “▶” to review changes on a tendency diagram with time, as well as changes of the tendency curves and values.
- (5) Press button “✕” to exit Tendency Diagram Review.

7.5Tendency Table Review

The tendency table data in the latest 120 hours can be displayed with resolutions as follows: 1min, 5min, 10min, 30min or 60min.

Directly click the option Review in the lower menu bar, or press the key “↶” to enter the Main Menu, enter the menu Review, and then select Tendency Table Review, as shown in the figure below:

TREND TABLE								
HR	60	60	60	60	60	60	60	60
SPO2	98	98	98	98	98	98	98	98
PR	60	60	60	60	60	60	60	60
NIBP	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)
IBP1	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)
	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)	110/80(95)
	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)	120/82(101)
	27	27	27	27	27	27	27	27
	39.0	39.0	39.0	39.0	39.0	39.0	39.0	39.0
	37.0	37.0	37.0	37.0	37.0	37.0	37.0	37.0
	09:31	09:32	09:33	09:34	09:35	09:36	09:37	
RES. [Icons]								

Figure 7-5 Tendency Table Review

See the relevant contents of “Tendency Chart Review” in the previous section to learn how to use the symbols in the interface.

(a) Tendency Table Review

Date is displayed below the measuring parameters and the time corresponding to each set of tendency

data is displayed on the right side of the displayed data. Parameters listed in a tendency table are divided into many categories, and the parameter of each instrument is displayed depending on the functional configuration of the instrument.

(b) Select Tendency Tables of Different Resolutions

Select Resolution to enter the menu and then select time interval for tendency data.

(c) Review Tendency Data of Different Parameters

Press “▲” “▼” to select and review other parameters not displayed on current views.

(d) Review Previous or Nearer Tendency Curves

Select “◀” “▶” to review previous or nearer tendency data.

(e) Operation Example

Review NIBP Tendency Chart:

- (1) Directly click the option Review in the lower menu bar, or press the key “↩” to enter the Main Menu, enter the menu Review, and then select Tendency Table Review;
- (2) Select the option Resolution in the window “Tendency Table Review”, and pick “1min”.
- (3) Select the key “▲” or “▼” (PageUp or PageDown) until NIBP tendency data is displayed on the screen.
- (4) Select “◀” “▶” to review previous or nearer tendency data.
- (5) Press the key “✕” to exit Tendency Table Review.

7.6 Waveform holographic Review

This monitor can display waveform holographic replay within 24 hours in the window Waveform Holographic Review (the specific time is related to the stored waveform and waveform quantity). You can review the waveform of any parameter for configured function of the monitor, as shown in the figure below:


Figure 7-6 Waveform Holographic Review

Please refer to Section 7.4 “Tendency Chart Review” to learn how to use the symbols.

(a) Operation Example

Observe the Waveform Review of ECG:

- (1) Directly click the option Review in the lower menu bar, or press the key “↶” to enter the Main Menu, enter the menu Review, and then select Waveform Holographic Review;
- (2) In the window Waveform Holographic Review, select the parameter to be reviewed, e.g. ECG waveform, select II or I in the Waveform Storage, and tick “√” in the front of the selected parameter.
- (3) In the window Waveform Holographic Review, press the key “▲” or “▼” (PageUp or PageDown) to select the desirable parameter;
- (4) Press “▶▶” to enter window on next page, and then press “◀◀” “▶▶” or “◀” “▶” buttons to review changes on a tendency diagram with time, as well as changes of the tendency curves and values.
- (5) Press the key “✕” to exit the window Waveform Holographic Review.

Chapter 8 Maintenance, Cleaning and Care

Use only materials and methods that are approved by our company and listed in this chapter for cleaning or disinfecting the device. Our company does not provide any guarantee for the damage caused by unauthorized materials or methods.

Our company has no responsibility for the effectiveness of controlling infectious diseases using these chemical agents. Please contact the infectious or epidemic disease experts in your hospital for details. Also see all policies which are suitable for your hospital and locality.

8.1 General

The monitor must be kept dust-free. After monitor cleaning and disinfection, please carefully check the monitor. If you find any signs of damage or wear on the monitor, stop using the monitor. If necessary, please clean it first before returning it to Shenzhen Comen Medical Instruments Co., Ltd. Please pay special attention to the following:

- Follow the manufacturer's instructions to dilute the solution, or adopt the lowest possible concentration.
- Do not let liquid enter the monitor.
- Do not pour liquid onto the monitor.
- No part of this monitor can be subjected to immersion in liquid.
- Don't use abrasive material (such as steel wool or silver polish etc) or bleaching powder, avoid using acetone-based cleaners such as acetone.

Warning

- Before cleaning the monitor or accessories, make sure that the equipment is switched off and disconnected from AC power.
- If any ECG cable is damaged or aging, the cable should be replaced with a new one.

Caution

- If you accidentally pour liquid onto the monitor or accessories, please contact our customer

service immediately.

 Attention

- Appropriate disinfection materials for the sensor, detector, cable, probe, are introduced in the product specification offered with the appendix.

8.2 Maintenance and check The overall check of the monitor, including a safety check, should be performed only by qualified personnel before first use, every 6 to 12 month, and each time after repair.

Before using the monitor, do the following:

- Check the work environment and if the power supply meets the requirement.
- Check if there is any mechanical damage.
- Check if the cables are worn and ensure insulation is in good condition.
- Check all the functions of the monitor to make sure that the monitor is in good condition.
- Check if the accessories used are specified by the manufacturers.
- Check the battery.
- If the monitor is equipped with a recorder, please check if the recorder is normal and recording paper meets the specified requirement.
- Check if the wiring resistance and leakage current meet the requirement.

If you find any damage on the monitor, stop using the monitor on patient, and contact the biomedical engineer of the hospital or our Customer Service immediately.

All the safety and maintenance checks that need to start the monitor should be performed by a qualified customer service technician. Non-professional operation can cause the monitor damage or cause a security risk, and human health may be endangered.

The circuit diagrams of the monitor can be provided by the manufacture, Shenzhen Comen Medical Instruments Co., Ltd as per customers demand. Qualified technicians can use it to help the user repair some apparatus that Shenzhen Comen Medical Instruments Co., Ltd classifies as “can be maintained by the user”.

 Warning

- If the hospital or agency that is responsible for using the monitor does not follow a satisfactory maintenance schedule, the monitor may become damaged, and human health may be endangered.

8.3 Maintenance schedule The following safety and maintenance check can be conducted by professional persons from our company. You can contact with our customer service technicians if you need the following maintenance check. Before the inspection or maintenance, the facilities should be cleaned and disinfected..

Check and maintenance	Frequency
According to IEC 60601-1 Medical electrical equipment Part 1 General requirement for safety	Check can be conducted at least every 2 years. Or after monitor falling,, power replacement, or when needed by customers
The ECG synchronism of the monitor and defibrillator	Check can be conducted at least every 2 years or when needed by customers.
NIBP air leakage check	Check can be conducted at least every 2 years or when needed by customers.
NIBP adjusting	Check can be conducted at least every 2 years or when needed by customers.
NIBP pressure calibration	Check can be conducted at least every 2 years or when needed by customers.
Mainstream and Side stream CO ₂ calibration and performance inspection	Check can be conducted at least every 2 years or when the measured value is in question.
Battery	See the section on battery for reference

8.4 Cleaning and Sterilization

Note on Cleaning and Sterilization

 Attention

- The monitor and accessory surface can be cleaned with hospital-grade ethanol and dried in air or with soft, clean cloth. For protecting the environment, disposable accessories should be recycled or disposed of properly.

 Caution

- Do not autoclave the sensor.
- Do not soak the sensor in liquids.
- If the sensor or cable is damaged or has signs of deterioration, do not use again.

8.4.2 Cleaning

The Monitor must be kept dust-free. Regular cleaning of the monitor shell and screen is strongly recommended. More cleaning is needed in the environmental polluted or sandstorm areas. Before cleaning the monitor or the sensor, please consult with customer service or understand the hospital equipment cleaning regulation.

(a) Cleaning agents are listed below:

- Diluted Ammonia Water
- Diluted Sodium Hypochlorite (Bleaching agent)
- Diluted Soapy Water
- Hydrogen Peroxide 3%
- Alcohol 70%
- Isopropanol 70%

(b) Before cleaning the monitor:

- a) Make sure that the equipment is switched off and disconnected from the power line.
- b) Use soft cotton ball to adsorb a small amount of cleaning agents and clean the screen.
- c) Use soft cloth to adsorb a small amount of cleaning agents and clean the monitor shell.

- d) If necessary, use a soft dry cloth to wipe away the excess cleaning agents.
- e) Dry the monitor in air.

8.4.3 Sterilization

To avoid damage to the monitor, sterilization is recommended only when stipulated as necessary in the Hospital Maintenance Schedule. Monitor facilities should be washed first.

Recommended sterilization material: Ethyl alcohol 70%, Isopropanol 70%, Glutaraldehyde 2%.

Caution

- Do not use ETO gas to disinfect the monitor.
- Use a moistened cloth to wipe up any agent remaining on the monitor.

Chapter 9 Patient Safety

9.1 Safety Instruction

The design of the patient Multiparameter monitor requires the international safety standard of medical electrical accessory. The instrument is provided by the anti-defibrillation with floating input and a protection of surgical electrotome. Use the correct electrode (refer to the chapter 'ECG Monitoring') and place it according to the guidance of the manufacturer. The display will resume 10 seconds after defibrillation.

9.2 Environment

The following guidance should be observed in the interest of absolute safety of electrical installations.

Vibration, dust, corrosives or explosive gas, extreme temperature and humidity should be avoided in the environment where the multi-parameter monitor is used.

Ventilation inside the instrument case where the monitor is installed should be guaranteed in such a manner that enough space is reserved in both the front to facilitate operations and the rear so that, when the case door is open, facilitate maintenance. A void space of at least 2 inches or 5 centimeters should be cleared around the instrument for air ventilation.

The monitoring system should be in the ambient temperature $-10^{\circ}\text{C} \sim +40^{\circ}\text{C}$ (keeping) $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$ (running) to satisfy the keeping and running requirements. An ambient environment out of this range may impair the instruments accuracy and cause damages to its components and circuits.

9.3 Power Requirement Please refer to the chapter of General Description or the chapter of Monitor Performance Indexes.

9.4 Grounding protection

To protect the patients and operation staff, the shell of patient multiparameter monitor must be grounded. So the monitor is equipped by a removable triaxial cable. When the cable plugs in a matching plug connector, the grounding wire of the power line ground the monitor. In case a three-plug connector is not

available, the electrical operating staff should be consulted.

 Warning

Replacement of the three-plug connector with a two-plug connector is strictly prohibited.

The grounding leads should be connected with the equipotential grounding terminal of the instrument. The instrument users do not know whether a given combination of instruments may invite dangers, e.g. due to accumulations of leaked currents, should consult relevant manufactures or experts in this field so as to guarantee that the required safety of the combined instruments are not compromised when the given combination is in use.

9.5 Equipotential Grounding

The primary protection of the instrument is embodied in the building protective grounding (protective ground) system by means of power plugs grounding. The multi-parameter monitor should be separately connected with the equipotential grounding system for examinations of hearts or skulls. One end of the equipotential grounding leads (potential equalizing leads) should be connected onto the equipotential grounding terminals on the rear panel of the instrument and the other should be connected onto one connector of the equipotential system. The equipotential grounding system should be in place for safety functions of the protective grounding leads in case of any damage to the protective grounding system. Cardiac or brain examinations should be conducted only in the rooms equipped with protective grounding systems. A check of the instruments should be conducted to guarantee the instruments are in good repair before each examination. The cables connecting with patients and instruments should be guaranteed not having been subjected to electrolytic pollution.

 Warning

- **Battery power should be used to power the monitor against unstable protective grounding (protective earthing) system.**

9.6 Condensation

The working instruments should be guaranteed not to form any condensation. Transferring of the instrument from one room to another may cause condensation on the instrument. This is attributed to its exposure to humid air at different temperatures. Unnecessary problems can be avoided by placing the

instrument in a dry place before putting it into use.

Note: Condensation is defined as coagulation of gases or liquids when cooled, e.g. water vapor when cooled is transformed into water and water when cooled into ice. The lower the temperature is, the faster condensation is formed.

Warning

- Use in the presence of combustible anesthetics is prohibited to avoid any risk of explosions.

9.7 Description of Symbols on Instrument Please refer to Section 1-1.6: Instrument and Symbols.

Chapter 10 ECG Monitoring

10.1 Definition of ECG monitoring

ECG produces a continuous record of electrical activity of the patients' heart and displays the activity on the monitor in the form of waves and values, so as to accurately assess the physiological status of the patient at that time. So, the connection of ECG cables should be assured to be accurate to obtain the correct measurements. In normal use, the instrument can only display one ECG waveform.

This ECG monitor works through C30 micro-host module. Regarding the monitors that is installed with micro-host, stick the micro-host into the monitor before starting monitoring.


Figure 010-11 micro-host cardial-electric cable portal

Patient cable consists of two parts:

- Trunk for connecting to monitor
- Lead device for connecting to patient

This monitor is applied in 3-lead, 5-lead and 12-lead monitoring. Different waveform displays will be acquired according to cardial-electric cables with different leads. Users can alter the lead names on the left of the ECG waveform on the screen to display the lead that requires emphasized monitoring. The monitor displaying parameters include heart rate (HR), ST segment analysis and ARR analysis.

The instrument has alarm functions for all of the above parameters.

 Attention

- Under factory settings in the instrument, the ECG waveforms are displayed in the top-two waveform positions.

10.2 Precautions for ECG Monitoring

 Warning

- ECG lead wire provided by the Comen Company should be used when monitoring the ECG signal by this instrument.
- When connecting the electrode or patient cable, you should insure that the patient is absolutely not in contact with any other conductive instrument or with the ground. In particular, you should make sure that all ECG electrodes, including the neutral electrodes, are connected to the patient, to prevent them contacting the ground or each other.
- Pacemaker failure: during the heart's complete conduction blockage or when Pacemaker can't wrestle/took away ,the high P waves (larger than the 1/5 of the average height of R wave) may be wrongly counted by the monitor, leading to asystoles being undetected
- Electrosurgery (ESU) equipment interference, defibrillation and discharge:
 - ✧ The operators shall not touch the patients, tables and instruments during the defibrillation period.
 - ✧ During the defibrillation, the electrocardio cable connected with patient may be damaged. Before reusing these cables, check whether the function is normal.
 - ✧ It can be recovered within 10 seconds after defibrillation, and will not lose any of the stored data. For electrosurgery (ESU) or defibrillation, measurement accuracy may temporarily drop, but it does not affect the safety of patients or equipment
 - ✧ Never expose the equipment under x-ray or strong magnetic field (MRI)

 Attention

- **Interference from ungrounded equipment near the patient and ESU interference may cause waveform problems. If you operate under the conditions regulated by EN60601-1-2 (with the anti-radiation capacity of 3V/m), the electric field strength more than 1V/m may cause measurement errors at various frequencies. Therefore, it is recommended not to use electrical radiation equipment near the ECG / respiration monitoring devices.**

10.3 Steps of Monitoring Preparation

Ask the patient to remove outer clothing on upper body:

- As skin is poor conductor of electricity, to get a good contact of the electrode and the skin, it is important to make preparation of the patient's skin.
- If necessary, shave the area for the electrode.
- Thoroughly clean the skin with soap and water. (Do not use ether or pure alcohol, because they will increase the resistance of the skin).
- Dry and rub the skin in order to increase the capillary blood flow and remove skin debris and oil.
- Attach the alligator clip prior to placement of the electrode.
- Place electrodes on the patient. If the electrodes used do not have conductive gel, apply the conductive gel before placement.
- Connect the electrode lead and the patient cable.
- Check that the monitor power is on.

Warning

- **ECG contact area should be checked daily for irritation. If there are signs of allergy, you should replace the electrode or change the position every 24 hours.**
- **Before ECG monitoring you should check whether the connections are normal or not. After pulling out the ECG cable, the screen will show "sensor fell off", and trigger the audible alarm at the same time**

 Attention

- In order to protect the environment, the electrodes used should be recycled or be disposed of properly.

10.3.2 Installing ECG lead

10.3.2.1 Identification and color code of electrodes

The table below shows the lead names in the European and American standards. (The leads are represented by RA, LA, RL, LL, and V in the American standard, while by R, L, N, F, and C in the European standard.):

Identification and color code of three-lead and five-lead electrodes are as follows:

American standard		European standard	
Name of lead	Color	Name of lead	Color
RA	White	R	Red
LA	Black	L	Yellow
LL	Red	F	Green
RL	Green	N	Black
V	brown	C	White

Identification and color code of twelve-lead electrodes are as follows:

lead	the American standard		the European standard	
	Marking	Color	Marking	Color
right arm	RA	white	R	red
left arm	LA	Black	L	yellow
right foot	RL	green	N or RF	Black
left foot	LL	red	(f)	green

chest 1	V1	red	C1	red
chest 2	V2	yellow	C2	yellow
chest 3	V3	green	C3	green
chest 4	V4	blue	C4	Brown
chest 5	V5	orange	C5	Black
chest 6	V6	purple	C6	purple

10.3.2.2 Installation position of three-lead ECG monitoring electrodesThe position of electrode of three lead device, is divided into American standard and European standardsee Figure10-2):

White / red (right arm) electrodes - placed under the clavicle, near the right shoulder.

Black / yellow (left arm) electrodes - placed under the clavicle, near the left shoulder.

Red / green (left leg) electrodes - placed in the left lower abdomen.


Figure 10-22 the position of electrodes of three lead device

10.3.2.3 Installation position of five-lead ECG monitoring electrodesThe position of electrode of five lead device, is divided into American standard and European standard (see Figure10-3):

White / red (right arm) electrodes - placed under the clavicle, near the right shoulder.

Black / yellow (left arm) electrodes - placed under the clavicle, near the left shoulder.

Green/ black (right leg) electrodes – placed at the right lower quadrant.

Red / green (left leg) electrodes - placed in the left lower abdomen.

Brown/white (chest) electrode-as Figure 10-4 on the chest wall


Figure 010-33 the position of electrode of five lead device

As with 5-lead devices, the chest lead electrodes should be put in one of the positions as follows (Figure10-4) :

V1: in the fourth intercostal space on the right edge of the sternum.

V2: in the fourth intercostal space on the left edge of the sternum.

V3: in the middle position of V2 and V4.

V4: in the fifth intercostal space of the left mid-clavicular line.

V5: at the left anterior axillary line, with horizontal position the same as the V4.

V6: at the left median axillary line, with horizontal position the same as the V4.

V3R-V6R: in the right chest wall, with its position corresponding to the left position.

VE: at the eminence of xiphoid.

V7: in the fifth intercostal space of left posterior axillary line on the back.

V7R: in the fifth intercostal space of right posterior axillary line on the back.


Figure 010-4 the position of Chest lead electrode of five lead device

10.3.2.4 Installation position of reformed twelve-lead ECG monitoring electrodes The position of electrode of twelve lead device after being improved, divided into American standard and European standard (see Figure 10-5/6).

White / red (right arm) electrodes - placed under the clavicle, near the right shoulder.

Black / yellow (left arm) electrodes - placed under the clavicle, near the left shoulder.

Green/ black (right leg) electrodes – placed at the right lower quadrant.

Red / green (left leg) electrodes - placed in the left lower abdomen.


Figure 010-5 Position of electrodes of twelve lead device

The chest electrode placement sites are generally six, taking spatium interosseum of costal bone as the locating basis, V1 ~ V6:

V1/C1: in the fourth intercostal space on the right edge of the sternum

V2/C2: in the fourth intercostal space on the left edge of the sternum

V3/C3: at the midpoint of the C2 and C4

V4/C4: at the intersection of the fifth intercostal space on the left edge of the sternum and the left mid-clavicular line

V5/C5: at the parallel previous axillary line of C4 level on the left

V6/C6: at the parallel median axillary line of C4 level on the left


Figure 010-6 Position of electrodes of twelve lead device

Attention

- In order to ensure the safety of patients, all leads must be connected to the patient.

10.3.2.5 ECG lead connection recommended for surgical patients**Warning**

- When using electricity surgical equipment (ESU), put the ECG electrode at the middle position between ESU earth plate and electricity surgical knife to avoid burn. the cable of electricity surgical equipment and the cable of ECG can't be entwined..
- In the use of surgical electric equipment (ESU), never let the electrodes be close to the electrosurgery equipment on the earth plate, or ECG signal would be interfered a lot.

The position of ECG lead depends on the type of surgery. For example, for the thoracotomy, electrode can be put in the chest or back side. in the operating room, because use electric knife equipment, sometimes knife surgery artifact may affect ECG waveform, in order to cut the artifact, you can put the electrode on left and right shoulder, near the abdomen of left or right side, and the breast lead can be put in the left side of chest to avoid the electrode on the upper arm, or ECG wave will become very small.

10.3.2.6 Waveform quality

The users can arrange the lead according to their needs. The name of lead on the channel is displayed on the left side of the corresponding waveform, you can directly select them and make changes. Two waveform

channels, only channel 1 can display the ECG waveform, click the screen ECG waveform , enter the ECG waveform settings, select the appropriate name of the lead from the “ I , II, III, the AVR, AVL, AVF, V” and select the way of gain and filtering according to your need.

Feature of a good signal

- Tall, narrow and without notch.
- R wave tall, completely above or below the baseline.
- T waves less than 1/3 height of the R wave.
- P wave should be much smaller than the T wave

In order to obtain 1 mV calibration ECG wave .ECG should be calibrated, at the same time, the screen displays the prompt: can not monitor the patient during calibration.


Figure 010-7 standard ECG wave form

Attention

- If the electrode paste is correct, while the ECG waveform is not accurate, then replace the lead.
- Interference of the instrument which is ungrounded and electrical surgical equipment (ESU) near the patient's may cause the waveform problems

10.4 ECG display and screen hotkey


Figure 11-8 ECG hotkey

Attention

- The pacing signal detected will show as “**!**” at the top of the waveform of the ECG waveform

- ① the name of ecgThe contents of the ECG 3-lead, 5-lead or 12-lead ECG, the ECG shows a different name, detailed for ECG waveform settings menu.
- ② Cardiac waves gain

You can select each calculation channel gain, including $\times 0.25$ $\times 0.5$ $\times 1$, $\times 2$ files and automatic mode. The automatic way means that the gain is automatically adjusted by the monitor. The left side of the each ECG waveform channel have given ruler 1mv. The height of the 1mv ruler is proportional to the amplitude

! Attention

- When the input signal is too large, the peak may be truncated. The user can refer to the actual waveform to manually change the gain of the ECG waveform in order to avoid the waveform display incomplete.

- ③ Filtering method:

The diagnostic mode is unfiltered ECG wave; the monitoring way may lead to the artifact filter out caused by false alarms; surgical method in the operating room can reduce the artifact and interference from electrical surgical equipment. Filtering methods simultaneously effect on two channels, and display at the top of the first ECG waveform.

! Warning

- ECG waveforms have varying degrees of distortion in the filter mode of “monitoring” and

“operation” . At this time the system can only provide the basic status of the ECG, and it may have greater impact for the result of the ST segment analysis. so we recommended that when the interference is small, try to use the diagnostic mode for the monitoring of the patient

④ 1 mV of the ruler stick

Cardiac waves gain.

⑤ the upper and lower limits of ECG alarm

Display the current setting of upper and lower limits of ECG alarm

⑥ PVCs (arrhythmia) and ST1/2 segment analysis

Display the current state of PVCs and ST1/2 analysis, and display according to 1 second refresh rate for each parameter in turn.

When the PVCs turn on, it can monitor arrhythmia.

When the ST1 / 2 segment analysis turn on, ST segment is under monitoring.

⑦ Data of ECG:

Shows the current measured data of ECG

⑧ Wave form of ECG:

Shows the current ECG waveform

10.5 Setting of ECG

10.5.1 Set the type of lead

This monitor is applicable to the 3-lead, 5-lead and 12-lead ECG cable. When ECG cables of different leads are used, the ECG that can be monitored is different.

When ECG uses the 3-lead, the monitoring lead include I, II, III.

When the ECG uses 5-lead, the monitoring lead include I, II, III, aVR, aVL, aVF, V.

When the ECG use 12 leads, monitoring lead include I, II, III, aVR, aVL, aVF, V1 and V2, and V3, V4, V5, and V6.

Set the instrument type which is equipped with lead of ECG cable:

- (a) In the conventional interface, click the ECG waveform ,in the ECG waveform pop-up menu, select [ECG settings]
- (b) Select the lead type] → [3 lead, [5 ECG] or [12-lead].
- (c) The setting is complete. Press “” key to exit.

10.5.2 Set the name of focus monitoring lead

In the conventional interface, using 3-lead monitoring can only display one ECG waveform; and using 5-lead or 12-lead ECG can show two ECG waves.

- (a) Select the first channel of the ECG waveform, in the ECG waveform pop-up menu, select [name of the lead, select the focus monitoring lead, such as conductivity [II], press “” key to exit.
- (b) If the 5-lead or 12-lead ECG, select the second channel of the ECG waveform, in the ECG waveform] pop-up menu, select [name of lead, select the focus monitoring lead, for example, [I]guide;
- (c) The setting is complete. Press “” key to exit.

10.5.3 Set the gain

Gain is used to adjust the size of ECG wave amplitude. You can select each calculation channel gain, gain include $\times 0.25$ $\times 0.5$ $\times 1$, $\times 2$ files and automatic mode, automatic way is automatically adjust the gain by the monitor. The right side of the each ECG waveform channel has given ruler 1mv. The height of the 1mv ruler is proportional to the amplitude.

- (a) Select an ECG waveform, in the ECG waveform pop-up menu, select [gain] → [$\times 0.25$], [$\times 0.5$], [$\times 1$], [$\times 2$] or [automatically]
- (b) The setting is complete. Press “” key to exit.

Attention

- When the input signal is too large, the peak may be truncated. The user can refer to the actual waveform to manually change the gain of the ECG waveform in order to avoid the waveform display incomplete.

10.5.4 Set the filtering method

There are three choices of filtering methods: The diagnostic mode is unfiltered ECG wave; the monitoring way may lead to the artifact filter out caused by false alarms. The surgical method in the operating room can reduce the artifact and interference from electrical surgical equipment. Filtering methods simultaneously effect on two channels, and display at the top of the first ECG waveform.

- (a) Select an ECG waveform, in the ECG waveform pop-up menu, select [filtering method] → [diagnosis] [monitoring 2] or [surgery];
- (b) The setting is complete. Press “” key to exit.

Warning

- Only in the diagnostic mode, System can provide the true signal which is untreated ECG waveforms have varying degrees of distortion in the filter mode of “monitoring” and “operation”. At this time the system can only provide the basic status of the ECG, and it may have greater impact on the result of the ST segment analysis. In the operation mode, the analysis result of AWRR also has some impact. So we recommended that when the interference is small, try to use the diagnostic mode to monitor the patient.

10.5.5 Set the calculation channel

Heart rate calculation channels:

[Channel 1] represents according to the first ECG wave waveform data to calculate the heart rate.

Channel 2] represents according to the second ECG wave waveform data to calculate the heart rate.

[Automatically] represents the monitor select to calculate the heart rate channel automatically.

- (a) Select the parameters of the [ECG], in the pop-up ECG settings] menu, select [calculation channel] → [Channel 1], [Channel 2] or [Automatic];
- (b) The setting is complete. Press “” key to exit

10.5.6 Set the power frequency inhibition

Power frequency inhibition, inhibit the signal collected 50Hz or 60Hz frequency componentsWhen the filtration method is not the diagnosis, the system automatically open the power frequency inhibition; filtering

method is the diagnosis, according to the need to select to open or close.

(a) Select the parameters of the [ECG], in the [ECG set] pop-up menu, set the [power frequency inhibition] as:

- [ON]: It is used when the waveform shakes become frequent (such as burr in waveform), and filters at the power supply frequency.
- [OFF]: No filtration.

(b) The setting is complete. Press "X" key to exit

10.5.7 ECG calibration

(a) Select the parameters of the [ECG], in the pop-up [ECG set] menu, select [ECG calibration], set the calibration status.

(b) To stop the calibration, return to the [ECG settings] menu and select [Stop ECG calibration].

When the ECG is in calibration, it can not monitor patients. At the bottom left of the screen, the instrument displays the prompt: can not monitor the patient during calibration

10.6 ST segment analysis

10.6.1 About ST segment analysis

The normal heart beat and atrioventricular pacing are used for ST segment analysis. The monitor analyzes these heart beats and calculates the elevation and depression of the ST segment. The monitor can display that information for ST numerical form. It can continuously monitor all available leads. For ST-segment analysis, the ECG waveform is not required to be displayed on the monitor. When the ST analysis is in the implementation, it usually uses a special filter which ensures diagnostic quality. If you use the "diagnosis" filter out of the ECG filtering mode to monitoring ECG, the ST segment of the ECG wave will be slightly different from ST segment in the ST fragment with the same waveform. In order to diagnostically assess the ST-segment, please always switch to filter mode of "diagnosis". You can also select the "monitoring" or "operation" mode, but ST-segment data will have a serious distortion.

ST segment analysis can measure the elevation or depression of the specified ST-segment of a lead.

Meaning of ST-segment measurement: A positive number indicates the elevation, a negative number

indicates depression

Measuring range of the ST segment: $-2.0 \sim +2.0$ mV,

10.6.2 The impact for ST-segment.

Some clinical conditions makes it difficult to get reliable ST monitoring, such as:

- Unable to obtain the low noise of the lead;
- the irregular baseline arrhythmia such as atrial fibrillation / atrial flutter;
- The patient is continuous ventricular pacing
- Patients with left bundle block

When these things happen, you should consider to close ST monitoring

Warning

- ST level change information provided by this monitor, its clinical significance should be decided by the physician

10.6.3 On/Off ST segment analysis

- (a) Select the [ECG] of parameter area, in the pop-up ECG set] menu, enter the menu [ST segment analysis];
- (b) Set the [ST segment analysis] to ON or OFF;
- (c) The setting is complete. Press “” key to exit

10.6.4 Point of ST segment analysis

Set the ST measurement point as R peak point, the ST measured values of every cardiac integrated wave, is the vertical distance difference between the peak point and the two measurement points, as shown below:

Figure 10-12 the ST analysis point

The ST measured values of every cardiac integrated wave, is the vertical distance difference between the peak point and the two measurement points, as Figure shown below:

Attention

- **If the patient's heart rate or ECG waveform changes obviously, you need to adjust the location of the point of ISO and ST. The abnormal QRS wave groups are not considered when the ST section analysis is made.**
- ❖ Method to adjust the ISO and ST point:
 - (a) Select the [ECG] of parameter area, in the pop-up ECG set] menu, enter the menu [ST segment analysis];
 - (b) Set the [ST segment analysis] to ON;
 - (c) Select the correct ST-segment analysis point to enter the [point of ST segment analysis] window;
 - (d) Select [ISO], [ST] point. Press ◀ and ▶ to make adjustment.
 - (e) The setting is complete. Press “Enter” key to exit

Warning

- **Please ensure that the location of the ST measurement points suitable for the patients.**

10.7 Arrhythmia analysis

Arrhythmia analysis is used for the clinical care of patients with ECG to detect changes in heart rate and ventricular premature beat, to save the arrhythmic events and to generate alarm information. Arrhythmia analysis can be used to monitor pacemaker and non-pacemaker patients. Doctors may, based on the arrhythmia analysis, evaluate the patient's condition (such as heart rate, the PVCS (ventricular premature

beat) frequency, rhythm, and abnormal heartbeat) and give the appropriate diagnosis and treatment. In addition to the detection of ECG changes, the arrhythmia analysis can also be used to monitor patients and provide appropriate warnings.

The arrhythmia monitoring is disabled by default for the monitor turned off. Users can turn on this function according to your need.

The function of arrhythmia monitoring can, by testing and classification of arrhythmias and heart abnormalities, remind doctors to pay attention to the rhythm of the patient, and then give alarm.

The monitor can make 13 kinds of arrhythmia analysis.

Arrhythmia analysis system will store 100 alarm events, the operator can edit the arrhythmic events in the menu.

10.7.1 On/Off arrhythmia analysis

- (a) Select the parameters of the [ECG], in the pop-up [ECG set] menu, enter the [arrhythmia analysis] menu;
- (b) Select [arrhythmia analysis], select [On] or [off];
- (c) The setting is complete. Press “✕” key to exit.

10.7.2 Arrhythmia alarm settings

- (a) Select the parameters of the [ECG], in the pop-up ECG settings] menu, enter [arrhythmia analysis] menu;
- (b) Select [ARR alarm settings], select [ALL] or separate arrhythmia [ASYSTOLE], [VFIB / VTAC], [the R ON T], [VT> 2], [COUPLET], [PVC], [BIGEMINY], [TRIGEMINY], [TACHY], [BRADY], [PNC], [PNP], [MISSED BEATS];
- (c) The setting is complete. Press “✕” key to exit.
- (d) Others:

Upper and lower limit of alarm:

Name	Maximum	Minimum	Single regulating variable
------	---------	---------	----------------------------

PVCs	10	1	1
------	----	---	---

10.7.3 Arrhythmia Review

- Select the parameters of the [ECG], in the pop-up [ECG set] menu, enter [arrhythmia analysis] menu;
- Select [arrhythmia review] to enter the review window, as shown below:


Figure 10-8 Arrhythmia review

- Select the [start time] and the type of event [ARR], usually, the event type is the acquiescent state, and can not be selected.
- Click the alarm events, such as “*** ventricular premature bigeminy” to enter another window to review the parameters and waveform of the alarm, as shown below:


Figure 010-9 Arrhythmia review

- (d) And you can select [the waveform holographic review] a comprehensive review of the waveform of the arrhythmia.

Chapter 11 RESP Monitor

11.1 RESP (Breathing) Measurement

11.1.1 Generation of Breath


The monitor measures breath according to the thoracic impedance values of the two electrodes. The impedance variations between the two electrodes (because of the thoracic activities) produce a breathing wave on the screen.

11.1.2 The Settings of Breathing Monitor

To monitor respiration, we don't need additional electrodes, but the position of the electrodes is quite important. For example:

■ Lateral thoracic expansion

Some patients' thorax may expand laterally, especially for the neonates. At the right time, two breathing electrodes shall be placed on the right mid-axillary line and the left lateral thorax where the patient has the maximum breathing movement to make sure that the breathing waveform is distinct

■ Abdominal respiration

For some patients the chest movement is limited so they primarily conduct a abdominal respiration And right now, you shall remove the electrode which was placed on the left leg to the left abdomen where there is maximum expansion to make sure the breathing waveform is clear.

⚠ Attention

- Breathing monitor can not be applied to the patient who moves frequently, because it may lead a wrong alarms.

11.1.3 Place Electrodes for the Breathing Measurement

For example, if the five leads are applied, the connecting methods are shown on the charts below (for the connecting methods of other leads, please refer to the related sections of Chapter 10):


Figure 11-1 Position of electrode of five lead device

⚠ Attention

- Put the green and the red electrodes at opposing angles so as to get the best breathing wave. You should avoid putting the electrodes over the liver area and the left ventricle on the line of the respiratory electrodes, which can help to avoid slight errors made from cardiac overlap or the pulsing blood flow. This is particularly important for newborn babies.

11.2 RESP Settings

11.2.1 Gain settings

Gain is used to adjust amplitude of breathing waveform. The alternative gains are $\times 1$, $\times 2$ gears.

- (a) Select RESP waveform, select[gain]→[$\times 1$]or[$\times 2$]in the shortcut[RESP waveform]menu.
- (b) Setting is complete. Press “**X**” to exit.

11.2.2 Set asphyxia warning time

- (a) Select PKW [RESP], in the shortcut [RESP setting] menu, select [asphyxia warning]→select [10s], [15s], [20s], [25s], [30s], [35s], [40s] or [none-warning]
 - [10s]: means asphyxia warning time, if the asphyxia time exceed the setting time and RESP [Alarmgrenze (off/on)] set as [on], the monitor will has a alarm and so on.
 - [none-alarm]: means asphyxia warning is shut down.
- (b) Setting is complete. Press “**X**” to exit.

11.3 Maintenance and Cleaning

Warning

- Before cleaning the monitor or the sensors, you must turn off the power and switch off the AC power.
- If any ECG cable is damaged or aging, the cable should be replaced with a new one.

Cleaning:

This monitor and the sensor surface can be cleaned with medical alcohol, air dried or dried with a clean and dry cloth.

Sterilization

In order to avoid long-term damage to the equipment, we advise you sterilize these products only when it is absolutely necessary in your instrument maintenance plan. We also suggest you first wash these sterilization instruments.

Recommended sterilization materials for monitor:

Ethanol: 70% alcohol, 70% isopropyl.

Acetaldehyde base

Sterilization

In order to avoid long-term damage to the products, we suggest you sterilize them only when your hospital regulations think it essential to do so. We also advise you first clean the sterilization products.

12.1 Definition of SpO₂ Monitoring

SpO₂ plethysmography parameters measure blood-oxygen saturation, that is, the percentage of the total oxyhemoglobin. When 97% of the total number of hemoglobin molecules combines with oxygen in the arterial blood's red blood cells, this blood will have 97% SpO₂ oxygen saturation, and at the same time the monitor reads 97% SpO₂ value. This value shows the percentage of oxygen-carrying hemoglobin molecules, constituting oxyhemoglobin. Furthermore parameters of the SpO₂ can also provide pulse rate signal and plethysmography wave.


Figure 12-1 SpO₂ probe port

12.1.1 Principle Of Measuring SpO₂ Plethysmography Parameter

Pulse oximetry is a measurement of oxygen saturation. It is a continuous, non-invasive way of determining the hemoglobin oxygenation saturation. It is involved in measuring how much light emitting from the sensor side passes through the patient's tissue (such as a finger or ear) and then reaches the other side of the receiver.

The sensor is usually able to gauge the 660nm wavelength of red LED, and 940nm of infrared LED While LED's maximum available output power is 4mW.

Though the amount of light passing through depends on many factors, most of them are constant. However, one of the factors is that blood flow in the arteries changes with time, since it is pulsatile. It is possible to

obtain arterial blood oxygen saturation by measuring the amount of light absorbed during the pulse. And surveying pulse itself can supply a “plethysmography” waveform and pulse rate signal.

The main screen can display “SpO₂” value and “plethysmography” waveform.

Warning

- If there is carbonyl hemoglobin, methemoglobin, or dye dilution chemical present, the SpO₂ value will deviate.

12.1.2 Recognition of the type of blood-oxygen

The types of oxygen in the arteries have been configured when the equipment is delivered from the factory, which is recognizable by the socket connected to the oxygen probe: Digital SPO₂: the blue socket, a circular single-slot one, is set for the blue oxygen sensor, whose color is the same as the socket. Nellcor SPO₂: the red socket, a circular single-slot one, is set for the red oxygen sensor, whose color is the same as the socket. Masimo SPO₂: the white socket, a circular single-slot one, is set for the white oxygen sensor, whose color is the same as the socket.

Warning

- This type of monitor can automatically identify the oxygen sensor. However, as the internal hardware has been fixed before delivery from factory, the use of improper oxygen sensor will lead to wrong outcomes of parameters.

12.1.3 Monitoring of Blood-Oxygen/ Pulse

Warning

- Do not twist the cable of electrical surgical equipment with the sensor's.
- Do not put a sensor on limbs with an arterial catheter or intravenous tube.

Attention

- Do not place the sensor on the body part doing measurements of blood pressure, because the process of measurement during blood pressure measurement may occlude the blood flow, and that will affect readings of oxygen saturation.

12.2 Attentions on Monitoring of SpO₂/ Pulse

Warning

- First check whether the sensor cable is normal before the monitor is started. When you unplug the cable of the SpO₂ sensor from the jack, the screen will display error message as “sensor off” and simultaneously trigger an audible alarm.
- If the sensor or its packaging has signs of damage, do not use it and return to the factory.
- Continuous and prolonged period of monitoring may increase the risk of undesirable changes in the skin characteristics, such as extremely sensitive, reddening, blistering, or even pressure necrosis etc, especially of the newborns, or patients with perfusion disorders and varying or immature forms of skin. Particular attention to checking placement of the sensor according to the changes in the quality of skin, correct optical alignment and attachment methods. Also check, periodically, the location where the sensor is attached and make a change of the position if there is a decline of skin quality. More frequent examination may be required due to different status of individual patients.
- To set the SpO₂ alarm limit as 100% is equal to turn off the upper limit alarm. High oxygen levels will have premature children in danger of infecting the crystal-like fibrous tissue disease. Therefore, the limit of oxygen saturation alarm must be chosen carefully based on generally accepted clinical practice.

Attention

- Make sure that the nails cover the light inside the probe. Its lines should be placed on the back of hands.
- SpO₂ value is always displayed in a fixed place.
- SpO₂ waveform and pulse volume are out of proportion.

12.3 Steps of Monitoring

Warning

- Select the appropriate placement according to the instrument and its supporting oxygen probe, which is fundamentally vital to a neonatal.

(1) SpO₂ plethysmography measurements:

- Turn on the monitor;
- Attach the sensor to the appropriate location of the patient's finger;
- Insert the connector of the sensor cable into the hole marked SpO₂ on its module.


Figure 0-4 Installation of the sensor

(2) Neonatal's SpO₂ plethysmography measurements:

The process of the neonatal's SpO₂ plethysmography measurement is basically the same as with adults. Following is an introduction to the neonatal oxygen probe and methods for its placement.

- (a) neonatal blood-oxygen probe Neonatal oxygen probe consists of the Y-shaped blood-oxygen probe and neonatal oxygen probe sheath. Insert the LED side of the Y-shaped probe in the upper groove of the sheath, and respectively the PD side of the probe within the lower (See Figure 12-2), then the neonatal blood-oxygen probe is shown in Figure 12-3.

Y-type blood oxygen probe

Blood oxygen probe sheath for newborns


Figure 0-5 Infant blood oxygen probe (1)


Figure 0-6 Infant blood oxygen probe (2)

- (b) Placement of the neonatal's blood-oxygen probe Clamp the probe to the newborn's hand or foot (as shown in Figure 12-4). Hold the sensor, pull the band and lay its 'V' shaped brim into the corresponding side of the 'V' shaped tank of the sheath, stretch the string to an appropriate length (20mm or so) and lay the other side's 'V' shaped edge of the rope into its corresponding the 'V' shaped groove, and then release the rope, with its both sides of the 'V' shape embedded in respective 'V' shaped groove, thread the first latch and lock it, shown in Figure 12-4. If the band is too long, thread the second latch. The necessity of locating the sensor in such a way has a purpose that the optical components will be in the correct and opposite position. Moreover, do not pull the band too tight, to prevent wrong measurements and seriously blockage of the blood circulation.


Figure 0-7 Installation of the infant blood oxygen probe

 Attention

- When the accurate positioning between the test site and the probe fails, it may result in wrong readings of blood-oxygen saturation, and even stop monitoring because of the failure of the search for the pulse wave. In this case you should re-position the two.
- Excessive movement of measured sites may affect the accuracy of the measurement, therefore, you should calm the patient or replace sites in order to reduce the impact of excessive movement.

 Warning

- In a long and continuous monitoring process, check the condition of the peripheral circulation and skin under measuring every 2 hours or so, and if negative conditions happen, timely change the site under measurement.
- In a long and continuous monitoring process, it is advisable to check periodically the positioning of the probe to avoid inaccurate measurement due to changing in the positioning from moving or other factors.

12.4 Limits in Measurement

During operation, the following factors can affect the accuracy of blood-oxygen saturation measurement :

- High-frequency radio interference, for instance, interference self-generated from the host system or from electrical scientific instruments connected to the system.

- During magnetic resonance imaging scanning (MRI), do not use the photoelectric oximeter and oxygen sensor, since induced currents may cause burning.
- Intravenous dye.
- Patient's excessive movement.
- External radiation.
- Improper installation of sensor or improper contact position with the object .
- Sensor's temperature (optimal temperature should be among 28 °C - 42 °C) .
- The sensor is placed on the limbs with a blood pressure cuff, arterial catheter, or the pipeline of body cavity.
- The concentrations of non-functional hemoglobin such as carboxyhemoglobin (COHb) and methemoglobin (MetHb) etc.
- Extremely low degree of blood-oxygen saturation.
- The measured area has poor circulation.
- Syndromes such as shock, a Extreme low degree of blood-oxygen saturation, anemia or low temperature etc and application of vasoconstrictor drugs can reduce blood flow to the level of not being able to be measured.
- Measurement also depends on both the oxyhemoglobin and reduced hemoglobin's absorption of specific wavelengths of light. If any other factors absorb the same wavelength, they will generate false measurement, lower SpO₂ values. These factors are as follows: carbonization of hemoglobin, methemoglobin, methylene blue , indigo rouge. It is recommended to use only the SpO₂ probe described in the accessories.

12.5 SpO₂ settings

12.5.1 Set intelligent alarm

- (a) Select [SpO₂] in the parameter zone, select [intelligent alarm] in the [SpO₂ setting] menu that popped up→select [10], [25], [50], [100]seconds or [Disable];
- (b) Upon completion of setting, press the key “✕” to exit.

Intelligent alarm is designed to reduce false alarms, help the doctor grasp more accurate and timely changes in blood-oxygen. For example, set the intelligent alarm range as 50, the upper limit of NELLCOR blood oxygen alarm as 97%, and lower limit as 90%, the measured blood oxygen value turns out to be 80% and remains for 3 seconds, then, it decreases to 78% and remains for 2 seconds, so calculating from the time the alarm limit is exceeded, the sound-light alarm will initialize instantly only when the blood oxygen value is out of the alarm range for 5 seconds. Meanwhile , the circle beside the blood oxygen value returns to origin. (This function is only effective for NELLCOR blood-oxygen) Its calculation is as follow:

$$\text{Minus percentage points} \times \text{seconds} = \text{SatSeconds integer}$$

values of the calculated SatSeconds are shown:

%SpO ₂	second	SatSeconds
(90%-80%) x	3 =	30
(90%-78%) x	2 =	24
Total SatSeconds	=	54


Figure 0-8 Example Figure

SatSeconds Example:

After about 4.9 seconds, the instrument will report SatSeconds alarm, because 54 is beyond the intelligent alarm range 50 SatSeconds.

Within a few seconds , saturation will fluctuate a bit, and not be stable. Generally, the patient's SpO₂ value

may fluctuate between the upper and lower limits of the alarm, and may re-enter to non-alarm range several times. In this volatile period, the system will store the positive and negative points of the SpO₂, till it reaches the SatSeconds limit or patient's SpO₂ values returned to and remains in the non-alarm range.

12.5.2 Set average time

The average time, SPO₂ value displayed on the monitor, is the result of averaged data collected during a period of time. The shorter the average time, then the faster is the monitors response to the patient's SPO₂ value changes, but lower measurement accuracy. Conversely, the longer it is, the more slowly the monitors response to the patient's SPO₂ value changes, but the higher the measurement accuracy will be. In the care of critically ill patients, to set a smaller average time is beneficial to the immediate analysis of the disease. This function is effective only with TMasimo blood oxygen) Setting procedure as follows:

- (a) Select [SpO₂]in the parameter zone, select [average time]in [SpO₂ setting]menu that popped up, →select [2s], [4s], [8s], [10s], [12s], [14s]or [16s];
- (b) Upon completion of setting, press the key “✕” to exit.

12.5.3 Set calculation sensitivity

Calculation sensitivity is classified as: normal, high and preventing the probe from falling off. According to the level, “avoid-probe off” has the highest sensitivity. For typical patient monitoring, use the “normal” sensitivity. As for those patients who have moist skin, active exercise , or for other reasons , the probe may be dropped off a patient's body, then use the “sensitive” sensitivity. If the patient has very low perfusion levels and wants to improve the sensitivity performance, please use the “avoid-probe off” sensitivity. This function is effective only with Masimo blood oxygen) Setting procedure as follows:

- (a) Select [SpO₂]in the parameter zone, select [calculation sensibility]in [SpO₂ setting]menu that popped up→select [normal], [high]or [preventing the probe from falling off];
- (b) Upon completion of setting, press the key “✕” to exit.

12.5.4 Set signal indication

Signal indication is referring to the case when the signal is turned on, there will be a signal collecting mark appearing under the SPO₂ waveform. The mark mainly indicates the quality of the signal in the collecting process. The signal disappears as it is turned off. (This function is effective only with Masimo blood oxygen) Setting procedure as follows:

- (a) Select [SpO₂] in the parameter zone, select [signal indication] in [SpO₂ setting] menu that popped up → select [on] or [off];
- (b) Upon completion of setting, press the key “” to exit.

12.6 Maintenance and Cleaning

Caution

- Do not autoclave the sensor.
- Do not soak the sensor in liquids.
- If the sensor or cable is damaged or has signs of deterioration, do not use it again.

Warning

- Before cleaning monitors or sensors, you must switch off and disconnect the AC power.

Cleaning:

Wipe the sensor's surface first with a cotton ball dipped alcohol or a soft cloth, then thoroughly dry cloth. The same method could be applied to the cleaning of the LED sensor and receiver devices.

Disinfect the cable with 3% hydrogen peroxide or 70 % isopropyl alcohol, and other active agents are equally effective, the joint can not be immersed in the above solution.

12.7 Masimo Information

■ Masimo Patent

It contains one or more of the following U.S. patents: RE38,492, RE38,476, 6,850, 787, 6,826,419, 6,816,741, 6,699,194, 6,684,090, 6,658,276, 6,654,624, 6,650,917, 6,643,530, 6,606,511, 6,584,336, 6,501,975, 6,463,311, 6,430,525, 6,360,114, 6,263,222, 6,236,872, 6,229,856, 6,206,830, 6,157,830, 6,067, 462, 6,011,986, 6,002,952, 5,919,134, 5,823,950, 5,769,785, 5,758,644, 5,685,299, 5,632,272,

5,490,505, 5,482,036, international patents or a item or a number of patents referred to in the www.masimo.com/patents. Including functions from products of Satshare ® and the U.S. Patent 6,770,028. Other patents are under application.

■ Other Information

© 2006 Masimo Corporation. Masimo, Radical, Discrete Saturation Transform, DST, Satshare, SET, LNOP, LNCS, and LNOPv are federal registered trademarks of Masimo Corporation's.

RadNet, Radicalsreen, signal IQ, FastSat, fastStart are trademarks of APOD and Masimo Corporation.

Chapter 13 NIBP Monitoring

In the adult and child modes, as compared with the results of inside-artery measurement or oscillation measurement conducted in typical patient groups, the results of blood pressure measurement determined by this equipment meet the requirements of “Dedicated Requirements for the Safety and Basic Performance of Automatically Circulated Noninvasive Blood Pressure Monitoring Equipment” (IEC 60601-2-30) in terms of average error and standard deviation. By referring to the auscultation method, the 5th Korotkoff tone is used to determine diastolic pressure.

In the neonate mode, as compared with the results of inside-artery measurement, the blood pressures measured by this equipment on the bodies of representative patients meet the requirements of “Dedicated Requirements for the Safety and Basic Performance of Automatically Circulated Noninvasive Blood Pressure Monitoring Equipment” (IEC 60601-2-30) in terms of average error and standard deviation.

13.1 General

Uses oscillation method for Non-invasive blood pressure (NIBP) measurement.

Available for adults, children and newborns. **When** measuring in child and newborn patients, you must ensure that the correct mode is selected (see menu settings of patient information). Wrong patient mode could jeopardize patient safety, because high blood pressure of adults does not apply to children and newborns.


Figure 13-1 NIBP Interface

13.2 NIBP Monitoring

13.2.1 NIBP measurement

Warning

- Before starting the measurements, make sure the selected criteria applies to your patients (adults, children, newborn).
- Do not install the cuff on a limb with an intravenous infusion tube or a catheter. During cuff inflation, if infusion is slowed down or clogged, the area around the tube may be damaged.
- Ensure that inflatable tubes that attached blood pressure cuff and monitor are smooth and with no tangles.
- No-invasive blood pressure measurement cannot be run in patients suffering from sickle-cell disease or skin damage or any time there may be anticipated damage.
- For a patients suffering from serious disturbance of blood coagulation mechanism, the decision to operate automatic blood pressure measurement must be made according to the clinical evaluation, because the friction area between body and sleeves has the risk of hematoma.

- (a) Insert the inflatable tubes in the blood pressure cuff on monitor interface, turn on instrument power supply.
- (b) Attach a blood pressure cuff to the patient's upper arm or thigh in the following ways, (Figure 13-1).
 - Verify that the cuff is completely deflated.
 - Use a size suitable cuff on the patient; make sure that the ϕ mark is placed right onto the appropriate artery. Be sure that the cuff is not too tightly wrapped around the limb, otherwise it may cause discoloration or even ischemia to the far end of the limb.


Figure 13-2 Cuff Usage

⚠ Attention

- The width of the cuff should be 40% (or 50% for neonates) of the circumference of the limb, or 2/3 of the length of the upper arm. The length of the inflated part of the cuff should be long enough to wrap 50~80% of the limb; a size unsuitable cuff will produce a wrong reading. If the size of the cuff has a problem, use a larger cuff in order to reduce errors.

Adult/neonatal/infant can use the cuffs repeatedly:

patient type	Limb Circumference	Cuff width	Inflatable pipe length
infant	10~19 cm	8 cm	2 m
children	18~26 cm	10.6 cm	
adult 1	25~35 cm	14 cm	
adult2	33~47 cm	17 cm	
leg	46~66 cm	21 cm	

Neonatal/infant/adult disposable cuffs:

Size	Limb Circumference	Cuff width	Inflatable pipe length
1	3.1~5.7 cm	2.5 cm	2 m
2	4.3~8.0 cm	3.2 cm	
3	5.8~10.9 cm	4.3 cm	
4	7.1~13.1 cm	5.1 cm	

- Check if the edges of the cuff are located inside the range marked by < - >. If not, replace the cuff with a more suitable one.

Connect the cuff and inflatable tube. Body parts used for pressure-measuring should be in the same horizontal location with patient's heart. If unable to do so, it is necessary to use the following correction method to modify the measurement results:

- If cuff is above the heart level location, 0.75mmHg (0.10kPa) should be added to the displayed value for per centimeter gap.
 - If cuff is below the heart level location, the displayed value should minus 0.75mmHg (0.10kPa) for per centimeter gap.
- (1) Confirm that the monitor method is correct (monitor method displays on the upper right corner of monitor information area), if you need to change the monitor method, please enter "NIBP setting" and change "measurement model".
 - (2) Press the " button in the front panel and start pressure inflation.

13.2.2 Action PromptOne automatic measurement

Enter into the [NIBP Settings] menu, check [Interval Time] ; the user can select the time interval for automatic measurement. Then press the " button in the Control Panel, the system will automatically conduct measurement in accordance with set interval time.

(2) Stop automatic measurement

At any point of the automatic measuring process, automatic measurement can be stopped by pressing the " button.

(3) Carry out a manual measurement

- Enter into the [NIBP Settings] menu, then press the " button on the control panel to start one-time manual measurement.
- After automatic measurement stops, press the " button and you'll start a manual measurement. If you press the " button again, the manual measurement will stop and continue with the automatic measurement.

(4) Conduct a manual measurement in automatic measurement process

Press the “” button on the control panel to conduct a manual measurement.

- (5) Stop manual measurement during measurement.

Press “” button on the control panel again, to stop a manual measurement.

- (6) For continuous measurement

Enter into the [NIBP Settings] menu, check [Continuous Measurement] to start continuous measurement. This process will last 5 minutes.

- (7) Pause in continuous measurement

At any point of the continuous measurement process, continuous measurement will be stopped by pressing the “” button.

Attention

- If you doubt reading accuracy, use the same method to check patient vital signs before checking the monitor function.

Warning

- If liquid splashes on the devices or accessories, especially when the liquid is likely to enter into channels or monitor, please contact the hospital's maintenance department.

13.2.3 Measurement Restrictions According to the patient's condition, oscillatory measurement has some restrictions. Such measurements are looking for regular impulse waves produced by arterial pressure. In the case the patients' condition makes this kind of detection difficult, measurement values become unreliable and load time increases. Users should be aware that the following conditions will interfere with the measurement method, so that pressure is not reliable or load time increases. In this case, the patient's condition will make measurement impossible.

(1) Patient Mobility

If patient is moving, shaking or in spasms, measurement will be unreliable even impossible, as these may interfere with the detection of the arterial pressure pulse and load time will be extended.

(2) Arrhythmia

If patient has shown arrhythmia caused by irregular heart beats, measurements are unreliable or even

impossible and load time will be extended.

(3) Heart-Lung Machine

If patient is connected to an artificial heart-lung machine, measure cannot be conducted.

(4) Pressure Change

If within a certain time, arterial pulse pressure is being analyzed to get the measurements, when blood pressure in patients is rapidly changing, measurement will be unreliable or even impossible.

(5) Severe Shock

If a patient is in serious shock or hypothermia, the pressure will be unreliable. Causes that reduce blood flowing to periphery would cause a decline in arterial pulse.

(6) Limit Heart Rate

Blood pressure measurement cannot be performed when heart rate is lower than 40bpm (beats per minute) or higher than 240bpm (beats per minute).

(7) Obese Patients

A thick layer of fat around a limb damps oscillations from the artery, thus preventing them from reaching the cuff. The accuracy is lower than the normal one.

13.2.4 NIBP Display

NIBP measurement results are displayed in the parameter area; the following figure is for reference only; the graphic displayed on your monitor might be somewhat different:

Figure 13-3 NIBP Measurement

1	Time of previous measurement	2	Measurement mode (patient type): Adult, child or neonate
3	Pressure unit: mmHg or kPa	4	Mean blood pressure
5	Cuff pressure displayed during measurement	6	Diastolic blood pressure
7	Measurement mode: Manual measurement, automatic measurement or continuous measurement	8	Systolic blood pressure
9	Interval time to the next automatic measurement	10	Displayed alarm limit of systolic pressure

13.3 NIBP Settings

13.3.1 Setting the Measurement Mode

Measurement modes include adult, child and neonate. Please select a mode suitable for the monitored patient. This mode is the same as the [Patient Type] set in [Patient Information].

- Select [NIBP] in the parameter area, in the [NIBP Settings] menu popping up, select [Measurement Mode] → select [Adult] , [Child] or [Neonate] ;
- Upon completion of setting, press the key “” to exit.

13.3.2 Setting the Interval Time

- Select [NIBP] in the parameter area, in the [NIBP Settings] menu popping up, select [Interval Time] → select [Manual] or [Automatic] [1min], [2min], [3min], [4min], [5min], [10min], [15min], [30min], [60min], [90min], [120min], [180min], [240min], [480min];

- Manual: Conduct one-time NIBP measurement manually under the manual mode when measurement is needed.

- Automatic: Under the automatic mode, the monitor conducts automatic, repeated NIBP measurement based on the set time interval.

(b) Upon completion of setting, press the key “✕” to exit.

13.3.3 Consecutive Measurement When the instrument is set to continuous measurement, it conducts automatic measurement continuously within 5 minutes.

Select [NIBP] in the parameter area, in the [NIBP Settings] menu popping up, select [Continuous Measurement], then the instrument starts measuring, and measures continuously within 5 minutes. If you want to end continuous measurement, you can press the blood pressure measurement button on the casing to stop this measurement.

Warning

- If under the automatic or continuous measurement mode the noninvasive blood pressure measurement time lasts too long, the limb rubbing the cuff may be accompanied by purpura, ischemia or neurotrosis. During patient monitoring, you should check the color, warmth and sensitivity of the far end of the limb from time to time. Once you observe any anomaly, you should place the cuff on another place or stop blood pressure measurement immediately.

13.3.4 Setting the Pressure Unit

- (a) Select [NIBP] in the parameter area, in the [NIBP Settings] menu popping up, select [Pressure Unit] → select [mmHg] or [KPa];
- (b) Upon completion of setting, press the key “✕” to exit.

13.3.5 Setting the Pre-Inflation Value

- (a) Select [NIBP] in the parameter area, in the [NIBP Settings] menu popping up, select [Pre-Inflation Value] → set a suitable pre-inflation value;
- Neonates: The range of the pre-inflation value is 60-120 mmHg, the default value is 100mmHg; Children: The range of the pre-inflation value is 80-200 mmHg, the default value is 120mmHg;
 - Adults: The range of the pre-inflation value is 80-240 mmHg, the default value is 160mmHg;

- (b) Upon completion of setting, press the key “” to exit.

13.4 NIBP Reset

Select [NIBP] in the parameter area, in the [NIBP Settings] menu popping up, simply select [Reset] . Reset can restore the inflation value of the blood pressure pump to the initial setting. If the blood pressure pump is working abnormally, it can be checked by means of reset, and will recover from any anomaly caused by an accidental reason.

13.5 NIBP Pressure Calibration

Manufacturers recommend a calibrated pressure gauge (or mercury sphygmomanometers) with a precision more than 1mmHg to be used for calibration. Select "calibration" to begin calibration and this item becomes "stop calibration", if you press the knob right now, the system will stop the calibration.

Warning

- **Calibration of NIBP measurements should be carried out once every two years (or by your hospital's maintenance of Constitution). Machine performance should be checked as following details .**

✧ Pressure sensor calibration steps:

Metal containers with a volume of 500ml $\pm$ 5% instead of cuffs. Connect a calibrated standard gauge with measurement error of less than 0.8mmHg, a spherical air pumps with a t-interface and inflatable tubes attached to NIBP jack on the module. Set the monitor to "calibration" mode, and then inflate the metal containers with ball type air pump to pressure of 0,50 and 200 mmHg respectively. The difference of pressure value between standard pressure gauge and monitor should be within the 3mmHg. Otherwise, please contact the Comen Company maintenance engineers.

Figure 13-4 NIBP Calibration Connection Schematic

13.6 NIBP Air Leakage Check

It is used to check if the NIBP measuring pump has a leak; when the NIBP cuff is put on, you can utilize this button to start the process of NIBP inflation, thus finding out if the sealing condition of the NIBP air circuit is good. If the leakage test is passed, the system will not give any prompt; if not, there will be a corresponding error prompt in the NIBP information area.

✧ Leak detection process:

- (1) Connect the cuff and NIBP to the opening of the monitor.
- (2) Wrap the cuff around a cylinder with proper size.
- (3) Enter “NIBP settings” menu.
- (4) Turn the knob, move the cursor to “leak testing”, press the knob. “Leak detection...” will display beneath the NIBP parameters area on the screen. The system has begun to perform leak detection.
- (5) system automatically inflated to the pressure of 180mmHg.
- (6) after about 20 seconds, the system will automatically open the valve, which means leakage measurement is complete.

If there is no prompt message in the NIBP parameter area, it means that the system has no sign of air leakage. If “Pump Leaking...” is displayed, it means that the air circuit may have an air leakage fault. Under such a circumstance, the operator should check the entire connection for loosening, and run the leakage test again after having made sure that the connection is errorless. If the fault prompt still appears, please contact the manufacturer for repairing.

Patient
Monitoring

Gas pipe

Figure 13-5 NIBP Air Leakage Checking Connection Schematic

 Warning

- This leak detection is different from those described in the standard EN 1060-1, which is for users to simply test air leakage in NIBP inflation. If the system displays the NIBP leaks at the end of testing, please contact the Comen Company maintenance engineers.

13.7 Maintenance and Cleaning

(1) Maintenance and Cleaning of Repeatedly Used Blood Pressure Cuffs

The cuffs can be placed in a conventional hot air oven for high pressure sterilization, disinfected by means of gas or radiation disinfection, or soaked in a decontamination solution for sterilization. But you should remember to remove the rubber bags when using such a method. The cuffs cannot be dry cleaned. The cuffs can be washed by machine or by hand; hand washing can prolong the service life. You should take out the latex rubber bags before washing. When the cuffs are completely dry after washing, you can put the rubber bags back in. You should reinstall the rubber bag into the cuff: First, place the rubber bag at the head end of the cuff, so that the rubber tube lines up with the big opening of the long end of the cuff, then curl up the rubber bag longitudinally and insert it into the big opening of the cuff, hold the leather tube and the cuff, shake the entire cuff until the rubber bag is set in place, insert the leather tube into the cuff and let it go out underneath the lining with small holes, as shown in the figure below:


Figure 13-6 Replacing the Rubber Band in the Cuff

(2) Disposable blood pressure cuffs

A disposable cuff can be used only for one patient. Do not use the same cuff for different patients. Disposable cuffs cannot be sterilized or disinfected by high pressure steam. Disposable cuffs can be washed with soap in order to control infection.

Chapter 14 TEMP Monitoring

14.1 Temperature Monitoring

This monitor has two temperature measurement channels; body temperature data can be measured by temperature probes.


Figure 14-1 Temperature Probe Interface

14.1.1 Temperature Measuring Setting

- (1) If you are using a disposable temperature probe, the temperature probe should be inserted into the socket, then connect the probe to the cable. For reusable temperature probe, it can be directly inserted into the socket.
- (2) Securely attach the temperature probe onto the patient's body.
- (3) Power up the system

Warning

- At the beginning of monitoring, you should check that the probe cable is normal. Then unplug the temperature probe cable from the jack, the screen will display the message "temperature sensor is off" and sound the alarm.
- Be careful in handling the temperature probe and cable, when not in use, the probe and cable

should be pulled into the loose ring. If the wire is pulled too tight, it will lead to mechanical damage.

- Calibrate the temperature measuring instrument at least once every two years (or according to the required time in the hospital directive rules). Please contact the manufacturer when calibration is needed.

Attention

- Disposable temperature probe can only be used once.
- During the monitoring process, the temperature measuring instrument will automatically check itself once per hour. Self-checking will last 2 seconds, and will not affect the normal working of the temperature monitor.

14.2 TEMP Unit Setting

- Select [TEMP] in the parameter area, in the [TEMP Settings] menu popping up, select [Temperature Unit] → select [°C] (Celsius degree) or [°F] (Fahrenheit degree) according to your habit;
- Upon completion of setting, press the key “✕” to exit.

14.3 TEMP Display

The parameter area of the monitor can display the values and units of two channels of body temperatures ([T1] and [T2]), as well as the difference ([TD]) of the two temperatures. By selecting [TEMP] in the parameter area, you can open the [TEMP Settings] menu.

14.4 Maintenance and Cleaning

Warning

- Turn off the power and disconnect the AC power supply before cleaning monitor or the connected sensor.

Maintenance and Cleaning of a Repeatedly Usable Temperature Probe:

- The temperature probe cannot be heated up to over 100°C (212°F). It can only tolerate a temperature somewhere between 80°C (176°F) and 100°C (212°F) for a short time.
- Probe cannot be disinfected by steam. It can only be disinfected by alcohol-containing detergent.
- When using straight beam probe, please use a protective cover to shield it.
- When cleaning the probe, one hand holds the head end, another hand with a damp lint free cloth scrubs the probe in the direction of the probe connector.

**Attention**

- If what you use is a disposable temperature probe, you are not allowed to sterilize or repeatedly use it.
- In order to protect the environment, disposable temperature probes should be reclaimed or appropriately disposed.

Chapter 15 CO₂ Monitoring

15.1 General

The contents of this section are just the descriptions of the ways of measurement by sidestream and mainstream CO₂ modules, which are different from the working mode of CO₂ measured in an anesthetic gas (AG). Please note their differences.

This instrument measures the CO₂ pressure of the patient's air circuit, can measure the content of end-of-respiration CO₂ (EtCO₂), the content of minimum CO₂ intake (Ins CO₂), and the airway respiration rate (AWRR), and display the CO₂ pressure waveform.

Monitors using the way of CO₂ measurement are divided into sidestream monitors and mainstream monitors.

- In the side-stream measurement mode, the respiratory gases through patient airways are sampled with a constant sampling flows and analyzed by remote CO₂ sensors built in the measurement system.
- In the mainstream measurement mode, CO₂ sensors are mounted on an airway joint that is directly inserted into the respiratory system of a patient.

CO₂ measurement provides:

- CO₂ waveform.
- End-of-respiration CO₂ (EtCO₂): CO₂ value measured at the end of respiration.
- Minimum CO₂ intake (Ins CO₂): The minimum value measured during respiration.
- Airway Respiration Rate (AWRR): respiration rate per minute derived from CO₂ waveform.

Warning

- The bumping and shaking of the CO₂ module should be avoided whenever possible.

Attention

- Do not use this instrument where combustible anesthetic gases are present in the environment.

- This instrument can be operated only by professionals who are occupationally trained and familiar with this manual.

In the “Monitor Settings” main menu, select “Module Switch Settings”, set the CO₂ switch to On, as shown in the figure below:

Figure 15-1 CO₂ Waveform and Parameter Display

15.2 Understanding CO₂ Modules


Figure 15-2 Minihost CO₂ Module Interface


Figure 15-3 CO₂ Plug-in Module

This monitor can measure by using several CO₂ modules designated by this company:

- By connecting to an external CO₂ analyzer, CO₂ measurements on the minihost can be directly used by this monitor.
- When this monitor has the CO₂ measurement function on the CO₂ module (as shown in Figure 15-3), the minihost and the AG module, it gives priority to the CO₂ module (as shown in Figure 15-3).
- If you use the AG module to measure CO₂, please read relevant contents in the AG section.

15.3 Measuring Principle and Working Process

The CO₂ measuring principle is mainly based on the characteristic that CO₂ can absorb the infrared rays having a wavelength of 4.3μm. The measuring method works as follows: Gaseous CO₂ is introduced to a measuring chamber of which one side is irradiated by infrared rays, and sensors are employed to measure the attenuation degrees of received infrared rays at the other side of the measuring chamber, and the attenuation degree is directly proportional to the CO₂ concentrations.

The comparison expression for the conversion between CO₂ partial pressure and CO₂ concentration is:

$$\text{CO}_2 \text{ Partial Pressure (mmHg)} = \text{CO}_2 \text{ Concentration (\%)} \times \text{Pamp (Ambient Pressure)}$$

For example: 5% CO₂ = 38mmHg at 760mmHg

5% CO₂ = 35mmHg at 700mmHg

CO₂ Module: adopting Autorun instruction measurement mode, and the waveform is sampled once in every 31 milliseconds.

15.4 Operating Instruction for CO₂ Connection

- (1) The schematic of connection of the mainstream module produced by the RESPIRONICS company is shown in the figure below:


Figure 15-4 Mainstream Mode CO₂ Connection

- (2) The schematic of connection of the sidestream module produced by the RESPIRONICS company is shown in the figure below:


Figure 15-5 Sidestream Mode CO₂ Connection Schematic

- (3) The schematic of connection of the ISA™ sidestream analyzer produced by the PHASEIN company is shown in the figure below:


Figure 15-6 ISA™ Sidestream Analyzer (ISA CO₂) CO₂ Connection Schematic

- (4) The schematic of connection of the IRMA™ mainstream analyzer produced by the PHASEIN company is shown in the figure below:


Figure 15-7 IRMA™ Mainstream Analyzer (IRMA CO₂) CO₂ Connection Schematic

The monitoring equipment produced by this company supports CO₂ measurement by using a sidestream or mainstream module produced by the IRONICES company, or an ISA™ sidestream analyzer (ISA CO₂ (CO₂) CAT. NO. 800101) produced by the PHASEIN company.

 Attention

- When an ISA™ sidestream analyzer (ISA CO₂ (CO₂) CAT. NO. 800101) produced by the PHASEIN company is used for monitoring CO₂, please refer to the contents of the section titled “17.6 Measuring Procedure and Before-Using Checking” for the measuring procedure, and change the procedure for setting the AG module menu to CO₂ menu setting.

 Warning

- Before use, please check airway joints. Do not use when visible damage or breaks are found on the airway adapter.
- When CO₂ is not used, it must be turned off, otherwise the CO₂ module will be in a working condition all the time.

15.5 CO₂ Measurement Procedures

Based on different CO₂ modules used, determine if it is necessary to set such menu items as “Oxygen Compensation”, “Balancing Gas”, “Altitude” and “Atmospheric Pressure” etc for the instrument. When an ISA™ sidestream analyzer (ISA CO₂ (CO₂) CAT. NO. 800101) produced by the PHASEIN company is used, if there are no such settings as altitude and atmospheric pressure etc, it means that such functions have been carried by the module itself, need not to be set manually. Please refer to the instruction manual coming with the module.

If you need the CO₂ alarm message function, you can set this function in “CO₂ Settings”.

15.6 Measuring Procedure of RESPIRONICS Branded Mainstream and

Sidestream Modules

The RESPIRONICS branded sidestream analyzer operating procedure is roughly the same as the mainstream analyzer operating procedure; please refer to the sidestream analyzer operating procedure for the mainstream analyzer operating procedure.

- (1) Start the host monitoring equipment (if a minihost is used, please start the minihost monitoring equipment at the same time).
- (2) Or insert the CO₂ plug-in module into the host monitoring equipment; the indicator of the CO₂ plug-in module will illuminate, which means that the module has been successfully connected to the host monitoring equipment; otherwise, please reinsert the CO₂ plug-in module.
- (3) Make connections according to the CO₂ module type or Figure 15-4 or Figure 15-5, and connect the CO₂ module interface cable to the CO₂ interface of the CO₂ plug-in module or the minihost.
- (4) Enter into the conventional screen of the host monitoring equipment, select [Exchange Waveform] to call out the “CO₂” waveform and parameters which you want to monitor, such as [CO₂] (this step can be skipped if the screen has already displayed the “CO₂” waveform and parameters).
- (5) When the CO₂ module is connected to the monitor, its module working mode is in the “Measurement” state; however, in order to make sure that it is in the correct working state, please do enter into the [CO₂ Settings] menu to set its “Working Mode” to the [Measurement] mode.
- (6) Set [CO₂ Switch] to [On].
- (7) In the [CO₂ Settings] menu, set [Oxygen Compensation] to 21 (usually, although it is 21 in this menu, in order to make sure that its datum is in an activated state, the customer still needs to reselect it);
- (8) In the [CO₂ Settings] menu, select an appropriate [Balancing Gas]: Indoor air, laughing gas or helium (usually if there is no unused laughing gas or helium indoors, you can just select indoor air);
- (9) In the [CO₂ Settings] menu, select a correct [Altitude]: 0~5029. 2m, instrument default: 0m; mainly refer to the following table for its standard:

Air Pressure Conversion Table – End-of-Respiration CO₂ Data Read By Basing the Standard on Altitude

Sea-Level Elevation		Atmospheric Pressure	5%CO ₂
Inch	m	mmHg	ETCO ₂ mmHg
Sea Level (0)	Sea Level (0)	760	38

500	152.4	745	37
750	228.6	738	37
1,000	304.8	731	37
1,500	457.2	717	36
2,000	609.6	704	35
2,500	762	690	35
3,000	914.9	677	34
3,500	1066.8	665	33
4,000	1219.2	652	33
4,500	1371.6	640	32
5,000	1524	628	31
5,500	1676.4	616	31
6,000	1828.8	604	30
6,500	1981.2	593	30
7,000	2133.6	581	29
7,500	2286	570	29
8,000	2438.4	560	28
8,500	2590.8	549	27
9,000	2743.2	539	27
10,000	3048	518	26
10,500	3200.4	509	25
11,000	3352.8	499	25
11,500	3505.2	490	24
12,000	3657.6	480	24
12,500	3810	471	24
13,000	3962.4	462	23
13,500	4114.8	454	23

14, 000	4267. 2	445	22
14, 500	4419. 6	437	22
15, 000	4572	428	21
15, 500	4724. 4	420	21
16, 000	4876. 8	412	21
16, 500	5029. 2	405	20
16, 800	5120. 6	400	20

Table 16-1

Note: It is assumed that the atmospheric pressure is 760mmHg and the ambient temperature is 0°C at the sea level. Calculation of Atmospheric Pressure: the sea-level based ambient temperature is assumed as 0°C. Refer to the above Table.

Warning

- By setting sea-level elevation, the monitor is not automatically changed with air pressure compensations. Correct sea-level elevation must be set before the first use of CO₂ Measurement Program. Improper setting of sea-level elevation will result in incorrect CO₂ readings. A 5% CO₂ deviation is generally generated corresponding to difference of each 1000m height.

(10) In the [CO₂ Settings] menu, select a correct [Atmospheric Pressure]: 405~760mmHg, the instrument default is 760 mmHg; when CO₂ value is on the high side or on the low side, select an appropriate atmospheric pressure based on the local condition by referring to the table above (since the atmospheric pressure and the altitude are corresponding, the atmospheric pressure can be adjusted only by setting the altitude.)

(11) In the [CO₂ Settings] menu, select [Zeroing]; after zeroing, the following prompt is displayed at the lower right corner of the screen: Zeroing...Please Wait for 30s; you can start measuring CO₂ only after the prompt disappears.

15.7 Measuring Procedure of PHASEIN Branded Sidestream and

Mainstream Analyzers

The PHASEIN branded sidestream analyzer operating procedure is roughly the same as the mainstream analyzer operating procedure; please refer to the sidestream analyzer operating procedure for the mainstream analyzer operating procedure.

15.7.1 Measurement Steps

If you want to set the host monitoring equipment in order to start gas analysis, please execute the following procedure:

- a) Start the host monitoring equipment (if a minihost is used, please start the minihost monitoring equipment at the same time).
- b) Or insert the CO₂ plug-in module into the host monitoring equipment; the indicator of the CO₂ plug-in module will illuminate, which means that the module has been successfully connected to the host monitoring equipment; otherwise, please reinsert the CO₂ plug-in module.
- c) Connect the Nomoline sampling tube to the input interface of the ISA analyzer (CO₂ module)
- d) Connect the interface cable of the ISA analyzer to the CO₂ interface of the CO₂ plug-in module or the minihost.
- e) Enter into the conventional screen of the host monitoring equipment, select [Exchange Waveform] to call out the “CO₂” waveform and parameters which you want to monitor, such as [CO₂] (this step can be skipped if the screen has already displayed the “CO₂” waveform and parameters).
- f) When the CO₂ module is connected to the monitor, its module working mode is in the “Measurement” state; however, in order to make sure that it is in the correct working state, please do enter into the [CO₂ Settings] menu to set its “Working Mode” to the [Measurement] mode.
- g) **Set [CO₂ Switch] to [On].**
- h) Set appropriate [Oxygen Compensation], [Laughing Gas Compensation].
- i) To connect the outlet of the sample gas to the discharge system, or to make the gas to flow back to the patient’s circuit.
- j) If it is green LED indication, ISA Analyzer is available.
- k) To carry out inspection before use according to the statement in the “Inspection Prior to Use (2)”.

- l) If the inspection is normal, start to monitor the Anaesthetic Gas.

15.7.2 Check before use Before connecting the Nomoline sampling pipe to the breathing circuit, carry out the following steps:

- a) Connect the sampling tube to the gas entrance interface (LEGI) of the ISA CO₂ module.
- b) Check whether the green light of LEGI is steadily on or not (The indication system is normal.) .
- c) Exhale to the sampling tube, check if a valid CO₂ waveform and value are displayed on the host monitoring equipment.
- d) Use the finger tip to block up the sampling pipe, and hold on for 10 minutes.
- e) Examine whether there is obstruction warning and if the LEGI shows a red flashing light.
- f) Under proper circumstances: Carry out enclosure check on the patient's circuit that is linked with the sampling pipe.

Attention

- If the interface shows the prompt of 'no connection to the oxygen sensor', please reinstall the oxygen sensor.
- The end of the gas circuit adapter which connects the gas sampling pipe should point upward so as to prevent the condensing water drops entering the gas sampling pipe and blocking it up.

Warning

- Hang the external CO₂ analyzer onto the CO₂ bracket on the rear casing of the instrument; prevent the dropping damage of the CO₂ module.
- Unless HME is used to protect the IRMA probe, the state indicating LED should face upward all the time during IRMA probe placement.
- Do not pull the cable of ISA By-flow Gas Analyzer.
- Do not operate the ISA By-flow Gas Analyzer in the environment beyond the designated working temperature.
- Make sure all connections are firm and reliable. Any leakage will result in the inclusion of

ambient air in the patients respiratory gas, which leads to a wrong reading.

15.8 CO₂ SettingTurning on/off CO₂

After the CO₂ module is connected to the monitor, you need to turn on CO₂ settings; when CO₂ measuring stops, in order to protect the working life of CO₂, please turn off CO₂ by executing the following procedure:

- Select [CO₂] in the parameter area, in the [CO₂ Settings] menu popping up, select [CO₂ Switch] → select [On] or [Off];
- Upon completion of setting, press the key “” to exit.

15.8.2 Setting the Pressure Unit

- Select [CO₂] in the parameter area, in the [CO₂ Settings] menu popping up, select [Pressure Unit] → select [mmHg] or [KPa];
- Upon completion of setting, press the key “” to exit.

15.8.3 Setting Gas Compensation

- Select [CO₂] in the parameter area, in the [CO₂ Settings] menu popping up, select [Oxygen Compensation] → select [High], [Medium] or [Low];
- Select [Laughing Gas Compensation] →select [On] or [Off] ;
- Upon completion of setting, press the key “” to exit.

Warning

- Please set the compensation of various gases based on actual circumstances; otherwise the measured result may seriously deviate from the actual value, which will result in a wrong diagnosis.

15.8.4 Setting the Suffocation Alarming Time

- Select [CO₂] in the parameter area, in the [CO₂ Settings] menu popping up, select [Suffocation Alarm]

→ select [10s], [15s], [20s], [25s], [30s], [35s], [40s] or [No Alarm];

- [10s]: It indicates the suffocation alarming time; when the time for the patient to have suffocation exceeds the set time, and the CO₂ [Alarming Limit (On/Off)] is set to [On], the monitor will trigger an alarm. and so on.

- [none-alarm]: means asphyxia warning is shut down.

(b) Upon completion of setting, press the key “” to exit.

15.9 Discharging Waste Gases

When nitrous oxide and/or an anesthetic gas is used, you should prevent these gases from polluting the operating room. Usually the gas discharging outlet should be connected to (via the gas discharging pipe connected to the sample gas outlet of the host equipment):

A discharging system (used for discharging collected gases) or the patient circuit (used for the back flowing of collected gases).

Warning

- **Anesthetics:** When an anesthetic which is being used or a patient who recently used an anesthetic is measured, the gas discharging hole on the module must be connected to a waste gas processing system or the patient circuit (on the anesthesia machine or the respirator), so as to prevent medical personnel from inhaling the anesthetic.

15.10 Maintenance and Cleaning of RESPIRONICS Branded Mainstream and Sdestream Modules

15.10.1 Common Cleaning Clean with cloth, optionally dipped with 70% isopropanol, aqueous solution containing 10% sodium hypochlorite (bleacher), sterilizing spray cleaner (such as Steris Coverage Spray HB), ammonia water or mild soap water. Before cleaning, wash the cloth with rinse water and then wring out and air-dry the washed cloth. Make sure that sensor windows are clean and are air-dried before being used repeatedly.

15.10.2 Airway Adapter for Cleaning Reusable Mainstream Sensor

Rinse it with warm soapy water first and soak it in liquid sterilizing fluid, such as 70% isopropanol, aqueous solution containing 10% sodium hypochlorite (bleacher), 2. 4% glutaraldehyde solvent, e. g. Cidex Plus or Steris System 1 or ammonia water. Wash it with clean water completely.

15.10.3 Method for Sterilizing Reusable Adapter

Autoclave Sterilizer-operable only to adult-used adapters.

Ethylene oxide (ETO)- Sterilizing for 1. 5h.

Soaking in Cidex Plus for 10h.

Soaking in Perasafe for 10h.

U. S. Steris System 1 pasteurizer.

Before adapters are reused, please make sure that the windows are dried without any residuals and that adapters stand intact during operation or cleaning/sterilization.

15.10.4 Sterilization Times for Reusable Airway Adapter

The reusable airway adapters can be reused 100 times if the above sterilization method is used.

15.10.5 Zeroing

Please zero before monitoring CO₂; zeroing is to eliminate the effect of baseline drifting on the results during measurement, thus ensuring the correctness of measured results.

Usually, the module will zero itself automatically when necessary. The user can zero the module manually when the user considers it necessary: Select [CO₂] in the parameter area, in the [CO₂ Settings] menu popping up, select [Zeroing] to zero the CO₂ module. During zeroing, make sure that the patient circuit is exposed to the ambient air (21% oxygen and 0% CO₂) for approximately 30 seconds; when the 30s zeroing prompt on

the screen ends, it means zeroing is completed.

15.11 PHASEIN Branded Mainstream and Sidestream Analyzer Related Information

15.11.1 Zeroing

An infrared gas analyzer needs to determine the zero reference level for CO₂ measurement. This zeroing standard is called as “zeroing” here.

Automatic Zeroing

The ISA sidestream gas analyzer execute zeroing automatically by switching the gas sample from the respiration circuit to the ambient air. To execute automatic zeroing once every 24 hours, the ISA sidestream gas analyzer takes less than 3 seconds. If the ISA sidestream gas analyzer is equipped with an oxygen sensor, automatic zeroing also includes the indoor air calibration of the oxygen sensor.

Manual Zeroing

Select [CO₂] in the parameter area, in the [CO₂ Settings] menu popping up, select [Zeroing] to zero the CO₂ module. During zeroing, make sure that the patient circuit is exposed to the ambient air (21% oxygen and 0% CO₂) for approximately 30 seconds; when this menu is in a non default (settable) condition, zeroing can be executed.

Warning

- Since successful zeroing requires that the gas analyzer exists in the ambient air (21% oxygen and 0% CO₂), you should make sure that the ISA is placed at a well ventilated position. Before and after executing the zeroing procedure, avoid breathing in the vicinity of the ISA sidestream gas analyzer.

15.11.2 CO₂ Module Lighting Information

Overview of States Indicated by LEGI:

Indicating Signal	State
Not blinking green light	The system is normal

Blinking green light	Zeroing
Not blinking red light	Sensor Error
Blinking red light	Checking the sampling tube

15.11.3Others

The operating methods for CO₂ analyzers and AG analyzers produced by the PHASEIN company are identical. Therefore, for such contents as “adverse effects on performance, safety warning messages, airway obstruction, consumables, safety symbol information, patent and trademark, analyzer maintenance and cleaning” with respect to the CO₂ analyzer described in this section, the user can refer to relevant contents in Chapter 16, Anesthetic Gas Measurement.

16.1 General

Anaesthetic Gas (AG) can be used to measure the Anaesthetic Gas and Respiratory Gas of the patient under anesthetization. The Anaesthetic Gasmodule provides the numerical values at the end of exhalation (Et) and the inhalation (Fi) of the following gases.

1. Carbon dioxide (CO_2) – The measured value is EtCO_2 (maximum exhaled value – maximum exhaled value detected during respiration).
2. O_2 – Oxygen.
3. Nitrous oxide (N_2O) – Laughing gas.
4. MAC – Lowest pulmonary alveolus concentration, which is the basic index indicating the depth of inhalation anesthesia.
5. Airway respiration rate (AWRR) – The times of respiration per minute.
6. End-of-respiration (Et) and inhaled (Fi) values of AA:

Halothane: HAL


Isoflurane: ISO

Enflurane: ENF

Sevoflurane: SEV

Desflurane: DES

16.2 Understanding the AG Module

Figure 16-1 AG and O₂ Module Interface

16.3 Measurement Principle of Anaesthetic Gas

The Respiratory Gas can be analyzed by various measurement principles. The Dispersive Infrared (DIR) Method or the Nondispersive Infrared (NDIR) Method are commonly used to insulate the absorption characteristics of the gas sample. The DIR Method is to use a single optical light filter and a prism or a diffraction grating to separate the wave length of each kind of anesthetic. And NDIR Method is to get the infrared light go through several narrow-band light filters and determine what kind of gas exists in the mixed gas.

The most commonly used gas analytical method is the one based on the medium of NDIR Method. The measurement principle is based on that many gases absorb the infrared energy of a certain wave length.

Usually the sidestream sampling infrared multi-gas analyzer continuously sucks in gas samples from the junction of the patient circuit and an airway device (e.g. a face mask, an airway tube or a throat mask tube). A modern analyzer's rate of gas sampling from the respiration circuit is somewhere between 50ml/min and 250ml/min. Via a small cup or a sample compartment, the sample gas passes through the infrared transmitter, the light filter and the infrared detector. The output signal of the infrared detector is in proportion to the infrared energy which has not been absorbed by the gas. In order to simultaneously quantify and identify multiple gases, such as laughing gas, carbon dioxide and five inhalational anesthetic gases, multiple light filters are needed. Detected signals are amplified, and converted by means of complex

algorithm executed by a microprocessor. You need to note that oxygen cannot be detected by mean of infrared spectrometry. The gas analyzer measures oxygen by applying auxiliary technology such as paramagnet or oxygen sensor.

16.4 Anesthetic Gas Display

In “Monitor Settings” in “Main Menu”, select “Module Switch Settings”, set “AG” to “On”, return to “Main Menu”, then select “Anesthetic Gas Screen” in “Working Screen”; refer to the figure below:

Figure 16-2 AG Monitoring

The AG module can display measured waveforms and parameter values on the screen of the monitor, which include:

- CO₂, O₂, N₂O and AA waveforms;
- AWRR: Airway respiration rate;
- MAC: Lowest effective pulmonary alveolus concentration;
- End-of-respiration (Et) and inhaled (Fi) values of O₂, O₂, N₂O and AA.

AA stands for one of the Anaesthetic Gases of Des (Desflurane), Iso (Isoflurane), Enf (Enflurane), Sev (Sevoflurane) and Hal (Hhalothane).

When the AG analyzer detects two or more anesthetic gases, the primary anesthetic gas is displayed in the waveform area of the screen, but the secondary anesthetic (s) is/are not displayed, only “Two anesthetic gases are detected” is prompted in the message display area.

 Attention

- Only the waveform and value of one anesthetic can be displayed at a time. When the second anesthetic gas appears, “Two anesthetic gases are detected” is prompted in the message display area.

16.5 MAC (Minimum Alveolar Concentration) Calculation

MAC (MAC: minimum alveolar concentration) is a standard for comparing the effect of the inspiratory Anaesthetic Gases. The MAC value stands for the density of Anaesthetic Gas in the pulmonary alveolus (one barometric pressure) which is the density that will not provoke muscular movement response for 50% people tested on standardized pain stimulation.

If the mechanism of determining MAC value is carried out in the host equipment, the algorithm applied in the calculation process must be recorded in detail. The following formula can be used to calculate and show the MAC value by the density of the (exhaled) air at the end of the exhalation.

$$MAC = \frac{\%Et(AA1)}{X(AA1)} + \frac{\%Et(AA2)}{X(AA2)} + \frac{\%Et(N2O)}{100}$$

X(AA): HAL=0.75%, ENF=1.7%, ISO=1.15%, SEV=2.05%, DES=6.0%

For example, when using one or several anesthetics, the anesthesia module measures and get that the patient's air at the end of the inhalation contains 4% DES, 0.5% HAL and 50% N₂O. So the MAC value is equal to

$$MAC = \frac{4.0\%}{6.0\%} + \frac{0.5\%}{0.75\%} + \frac{50\%}{100\%} = 1.84$$

 Attention

- The above formula does not take the altitude, the patient's age and other personal factors into consideration.

16.6 Oxygen Sensors

● Paramagnetic Oxygen Analysis

The principle of paramagnetic oxygen analysis is to measure the attractive force produced by the high

intensity magnetic field to oxygen molecules in the mixed gas. A paramagnetic analyzer distinguishes oxygen from other gases based on the magnetic field magnetization sensitivity. Since oxygen has a paramagnetic property, it will be attracted by a magnetic field, while most of other gases will not be attracted by a magnetic field. In term of percentage, if oxygen is designated to a value of 100, the magnetic field magnetization sensitivity of most of other gases is close to zero.

Oxygen sensors applicable to ISA gas analyzers are PM1116 oxygen sensors made by Servomex. Such a sensor will produce a symmetric non-uniform magnetic field. If oxygen exists, it will be attracted to the strongest part of the magnetic force in that magnetic field. In that magnetic field, two glass balls filled with nitrogen are installed on a rotary suspension rod. A mirror is installed in the middle of the suspension rod. Light beams projected to this mirror are reflected to a pair of phototubes. The oxygen attracted to the magnetic field will push the glass balls from the strongest part of the magnetic force in the magnetic field, so that the suspension rod starts rotating. When the phototubes detect such rotation, they will generate a corresponding signal and transmit it to a feedback system. The feedback system will let an electric current pass through the cable on the suspension rod to produce a counteracting moment which keeps the suspension rod at its original position. This electric current is in proportion to the oxygen concentration, thus the oxygen concentration can be known by measuring the electric current passing through the cable.

The main advantages of the Paramagnetic Oxygen Sensors include:

- short rise time
- stability and precision
- No need to change or supply chemicals
- Generally dispense with maintenance

16.7 Measurement Steps and Inspection Prior to Use

(i) Measurement Steps

If you want to set the host monitoring equipment in order to start gas analysis, please execute the following procedure:

- (a) Insert the AG module into the host monitoring equipment.
- (b) To connect the Nomoline sampling pipe to the input interface of the ISA (the anesthesia module) Analyzer.

- (c) Connect the interface cable of the ISA analyzer to the host monitoring equipment.
- (d) Start the host monitoring equipment.
- (e) Enter into the [Measurement Settings] window of the host monitoring equipment, check if the "AG Module" is successfully connected; if the host monitoring equipment can display the AG module, it means that the AG module has been successfully connected.
- (f) Enter into the conventional screen of the host monitoring equipment, select [Exchange Waveform] to call out the "Anesthetic Gas" waveform and parameters which you want to monitor, such as [AA], [O₂], [N₂O], [CO₂ (AG)].
- (g) When the anesthesia module is connected to the monitor, its module working mode is in the "Measurement" state; however, in order to make sure that it is in the correct working state, please do enter into the [AA Settings], [O₂ Settings], [N₂O Settings] or [CO₂ Settings] menu to set its "Working Mode" to the [Measurement] mode.
- (h) To connect the outlet of the sample gas to the discharge system, or to make the gas to flow back to the patient's circuit.
- (i) If it is green LED indication, ISA Analyzer is available.
- (j) To carry out inspection before use according to the statement in the "Inspection Prior to Use (2)".
- (k) If the inspection is normal, start to monitor the Anaesthetic Gas.



Figure 16-3 ISATM Sidestream Analyzer Connection Schematic



Figure 16-4 IRMA™ mainstream Analyzer Connection Schematic


(j) Check before use

Before connecting the Nomoline sampling pipe to the breathing circuit, carry out the following steps:

- Connect the sampling pipe to the interface of gas inlet (LEGI) of the ISA Anesthesia Module.
- Check whether the green light of LEGI is steadily on or not (The indication system is normal.).
- As for the ISA Anesthesia Module furnished with oxygen option: Check whether the oxygen reading on the host equipment is correct or not (21%).
- Exhale to the sampling tube, check if a valid CO₂ waveform and value are displayed on the host monitoring equipment.
- Use the finger tip to block up the sampling pipe, and hold on for 10 minutes.
- Examine whether there is obstruction warning and if the LEGI shows a red flashing light.
- Under proper circumstances: Carry out enclosure check on the patient's circuit that is linked with the sampling pipe.

 Attention

- If the interface shows the prompt of 'no connection to the oxygen sensor', please reinstall the oxygen sensor.
- The end of the gas circuit adapter which connects the gas sampling pipe should point upward so as to prevent the condensing water drops entering the gas sampling pipe and blocking it up.

**Warning**

- Hang the external AG analyzer onto the AG bracket on the rear casing of the instrument; prevent the dropping damage of the anesthesia module.
- Unless HME is used to protect the IRMA probe, the state indicating LED should face upward all the time during IRMA probe placement.
- Do not pull the cable of ISA By-flow Gas Analyzer.
- Do not operate the ISA By-flow Gas Analyzer in the environment beyond the designated working temperature.
- Make sure all connections are firm and reliable. Any leakage will result in the inclusion of ambient air in the patients respiratory gas, which leads to a wrong reading.

16.8 Gas Settings

In the lower menu column or the “Main Menu”, select the “Measurement Settings” menu, enter into “Anesthetic Gas Settings”; you can set “Working Mode, CO₂, O₂, N₂O, Anesthetic Gases and Anesthetic Gas Zeroing”.

16.9 Setting AG

Select [AA], [O₂], [N₂O] and [CO₂ (AG)] in the parameter area; the [AA Settings], [O₂ Settings], [N₂O Settings] and [CO₂ Settings] menus pop up respectively.

16.9.1 Setting the Anesthetic Gas Working Mode

- (a) Select [AA] in the parameter area, in the [AA Settings] menu popping up, select [Working Mode] → [Standby] or [Measurement].
- Standby: When the anesthesia module is in a standby condition, the air pump is in a stop-working condition in order to prolong the module’s life; and “The anesthetic module is standing by” is prompted in the message display area.
 - Measurement: When the anesthesia module is in a working condition, the green light of the anesthesia module illuminates without blinking, and the anesthesia module sends measurement data to the monitor.

Normally, when the user connects the anesthesia module to the monitor on the screen, the module will automatically change its working mode to the “Measurement” condition, but the user must make sure that the anesthesia module's working mode is in a measurement condition.

- (b) Upon completion of setting, press the key “ $\times$ ” to exit.

16.9.2 Setting Gas Compensation

This option is in default state because the oxygen sensor provides oxygen compensation automatically according to the altitude. When anaerobic module provides oxygen compensation, compensation level can be set manually according to the altitude

- (a) For example, Select [AA] in the parameter area, in the [AA Settings] menu popping up, select [Oxygen Compensation] → [High], [Medium] or [Low].

- ✧ High: The system default oxygen compensation is 85%;
- ✧ Medium: The system default oxygen compensation is 50%;
- ✧ Low: The system default oxygen compensation is 21%;

- (b) For example, Select [CO₂ (AG)] in the parameter area, in the [CO₂ Settings] menu popping up, select [Laughing Gas Compensation] → select [On] or [Off].

- Laughing gas compensation: It is compensated automatically by the anesthesia module at starting.

- (c) Upon completion of setting, press the key “ $\times$ ” to exit.

Warning

- Please set the compensation of various gases based on actual circumstances; otherwise the measured result may seriously deviate from the actual value, which will result in a wrong diagnosis.

16.9.3 Zeroing.

An infrared gas analyzer needs to determine the zero reference level for CO₂, N₂O and anesthetic gas measurement. This zeroing standard is called as “zeroing” here.

- Automatic Zeroing

The ISA sidestream gas analyzer execute zeroing automatically by switching the gas sample from the respiration circuit to the ambient air. To execute automatic zeroing once every 24 hours, the ISA CO₂ gas analyzer takes less than 3 seconds, while ISA multi-gas analyzer takes less than 10 seconds. If the ISA sidestream gas analyzer is equipped with an oxygen sensor, automatic zeroing also includes the indoor air calibration of the oxygen sensor.

- Manual Zeroing

Select [AA] in the parameter area, in the [AA Settings] menu popping up, select [Zeroing]. During zeroing, make sure that the patient circuit is exposed to the ambient air (21% oxygen and 0% CO₂) for approximately 30 seconds; when this menu is in a non default (settable) condition, zeroing can be executed.

 Warning

- Since successful zeroing requires that the gas analyzer exists in the ambient air (21% oxygen and 0% CO₂), you should make sure that the ISA is placed at a well ventilated position. Before and after executing the zeroing procedure, avoid breathing in the vicinity of the ISA sidestream gas analyzer.

16.9.4 Setting the Pressure Unit

- Select [CO₂ (AG)] in the parameter area, in the [CO₂ Settings] menu popping up, select [Pressure Unit] → select [mmHg] or [KPa];
- Upon completion of setting, press the key “” to exit.

16.9.5 Setting the Suffocation Alarming Time

- Select [CO₂ (AG)] in the parameter area, in the [CO₂ Settings] menu popping up, select [Suffocation Alarm] → select [10s], [15s], [20s], [25s], [30s], [35s], [40s] or [No Alarm];
- [10s]: It indicates the suffocation alarming time; when the time for the patient to have suffocation exceeds the set time, and the CO₂ [Alarming Limit (On/Off)] is set to [On], the monitor will trigger an alarm. and so on.
- [none-alarm]: means asphyxia warning is shut down.
- Upon completion of setting, press the key “” to exit.

16.10 Anesthesia Module Lighting Information

Overview of States Indicated by LEG1:

Indicating Signal	State
Not blinking green light	The system is normal
Blinking green light ¹⁾	Zeroing
Not blinking blue light ²⁾	There is an anesthetic gas
Not blinking red light	Sensor Error
Blinking red light	Checking the sampling tube

Note 1: Not applicable to IRMA OR.

Note 2: Only applicable to IRMA multi-gas probes.

16.11 Adverse Effects on Performance

(1) The following circumstances can cause known adverse effects on indicated performance:

- Quantitative effects of humidity or condensation.
- Quantitative effects of atmospheric pressure.
- Disturbing gases or water vapor.
- Other disturbing sources.

(2) Gas measurement unit

Volume percentage is used as the unit of concentration of a reported gas. Concentration is defined below:

$$\%gas = \frac{\text{Partial pressure of gas component}}{\text{Total pressure of gas mixture}} \cdot 100$$

The total pressure of a mixed gas is measured by the cup pressure sensor in the ISA gas analyzer.

If you want to convert it into another unit, you can the actual atmospheric pressure sent by the ISA sidestream analyzer; for example:

$\text{CO}_2 \text{ (mmHg)} = (\text{CO}_2) \times (\text{atmospheric pressure (kPa) from the ISA}) \times (750 / 100)$.

For example: 5.0 vol% CO_2 @ 101.3 kPa $0.05 \times 101.3 \times 750 / 100 = 38 \text{ mmHg}$

(3) Effects of Moisture

The partial pressure and volume percentage of CO_2 , N_2O , oxygen or an anesthetic gas depend on the content of water vapor in the measured gas. Oxygen measurement will be calibrated; at the actual ambient temperature and humidity level, 20.8 vol%, instead of the actual partial pressure, will be displayed. 20.8 vol% of oxygen corresponds to the actual oxygen concentration of indoor air (the water concentration is 0.7 vol%) (for example: at 1013 hPa, corresponding to 25°C and 23% RH). The actual partial pressure at the current humidity level will be displayed all the time when CO_2 , N_2O , oxygen and anesthetic gases (all the gases measured the infrared cell) are measured.

In the pulmonary alveoli of a patient, the water vapor in breathing gases become saturated at the body temperature.

After breathing gases is collected and sent into the sampling tube, their temperature gets near to the ambient temperature before they enter into the ISA sidestream gas analyzer. After Nomoline has removed all the condensed water, water molecules will not enter into the ISA sidestream gas analyzer. The relative humidity of all the gases collected is approximately 95%.

If you need the CO_2 value under BTPS, you can use the following equation:

$$Et\text{CO}_2(\text{BTPS}) = Et\text{CO}_2 \cdot \left(1 - \left(\frac{3.8}{P_{\text{amb}}} \right) \right)$$

Where:

$Et \text{ CO}_2$ = $Et \text{ CO}_2$ value [vol %] sent by the ISA

P_{amb} = Atmospheric pressure [kPa] sent by the ISA

3.8 = Typical partial pressure [kPa] of water vapor condensed between the patient's circuit and the ISA

$Et \text{ CO}_2 \text{ (BTPS)}$ = $Et \text{ CO}_2$ gas concentration [vol%] under BTPS

It is assumed that standard oxygen of indoor air has been used at a humidity level of 0.7 vol% H_2O

16.12 Safety Alarm Information

16.12.1 ISA Sidestream Gas Analyzer Safety Warning Information

Warning

- The ISA sidestream gas analyzer is designed to be used by authorized or trained medical personnel.
- Only Nomoline sampling tubes produced by PHASEIN can be used.
- The ISA sidestream gas analyzer shall not be used in an inflammable anesthetic gas.
- You should earnestly neaten the sampling tube in order to reduce the risk of it wrapping or reining the patient.
- Do not repeatedly use a disposable sampling tube.
- Do not lift the ISA/host equipment by grasping the sampling tube; otherwise it may break away from the ISA/host equipment, which may result in that the ISA/host equipment falls onto the patient.
- Used disposable sampling tubes should be disposed according to local medical waste stipulations.
- Do not use a sampling tube configured for an adult or a child on an infant; otherwise the invalid cavity in the patient circuit will increase.
- Do not use a sampling tube configured for an infant on an adult; otherwise it will result in a too big flowing resistance.
- Do not use the ISA sidestream gas analyzer together with a quantitative spraying agent or spray; otherwise it may result in the clogging of the germ filter.
- Check if the flowing speed of the gas sample is too high for the given patient type.
- Since successful zeroing requires that the gas analyzer exists in the ambient air (21% oxygen and 0% CO₂), you should make sure that the ISA is placed at a well ventilated position. Before and after executing the zeroing procedure, avoid breathing in the vicinity of the ISA sidestream gas analyzer.
- The Nomoline sampling tube and its interface are not germ free devices. In order to prevent the sampling tube from causing damages, please never carry out high pressure sterilization on any part of the sampling tube.
- Never disinfect the ISA sidestream gas analyzer or soak it into a liquid.

- Mobile and radio frequency communication equipment will affect measurement. Make sure that the ISA gas analyzer is used in the electromagnetic environment designated in this operating instruction manual.
- The ISA gas analyzer can only be used as a piece of auxiliary equipment for patient evaluation. It must be used together with other vital sign and symptom evaluation equipment.
- If the input interface of the sampling tube starts showing red blinking, or a Nomoline clogging message is displayed on the host, the sampling tube should be replaced.
- It is not allowed to alter this equipment without the manufacturer's authorization. If this equipment is altered, appropriate checking and testing must be conducted in order to make sure that it can be safely operated over a long period of time.
- The ISA gas analyzer is not designed for being used in a MRI environment.
- During MRI scanning, the host equipment must be placed outside the MRI room.
- Using high frequency electrosurgical equipment in the vicinity of the ISA/host equipment may produce interference, which will result into incorrect measurements.
- Do not use the external natural heat dissipation function of the ISA equipment.
- Do not apply a negative pressure (such as using a syringe) onto the Nomoline to remove condensed water.
- If the positive or negative pressure in the patient circuit is too high, it may affect the sample's flowing speed.
- If the discharging or sucking pressure is too high, it may affect the sample's flowing speed.
- The discharged gas should be discharged into the patient circuit, or into a discharging system.
- If the collected gas sample needs to supply air for respiration, a germ filter should be used at the discharging side all the time.
- When placing the ISA gas analyzer, try not to place it at a position where the analyzer might fall onto the patient's body.

16.12.2 ISA Mainstream Gas Analyzer Safety Warning Information

Warning

- Do not use an IRMA airway adapter configured for an adult or a child on an infant, since the

adapter will add 6ml of invalid cavity to the patient circuit.

- If the airway adapter has water drops/condensation, it should be replaced.
- Use an IRMA airway adapter made by PHASEIN.
- Do not use an IRMA infant airway adapter on an adult; otherwise it will result in a too big flowing resistance.
- When an energized part is contacted, sufficient protection should be provided to the host equipment.
- Only an adapter cable approved by PHASEIN AB can be used.
- A warning must be implemented in the host equipment, displayed during demonstrative data displaying.
- The host equipment should be equipped with an appropriate alarm system to remind the user of circumstances which may cause death or serious damages to the patient's health.
- Every corresponding alarm message in IRMA state abstract fields must be implemented in the host equipment.
- The IRMA probe is not designed to be contactable to the patient.
- Incorrect probe zeroing will result in false gas readings.
- The IRMA probe is designed for being used by authorized or trained medical personnel.
- The IRMA probe is not designed shall not be used in an inflammable anesthetic gas.
- A disposable IRMA airway adapter shall not be used repeatedly. Repeatedly using a disposable adapter will cause cross infection.
- Used disposable airway adapters should be disposed according to local medical waste stipulations.
- Only oxygen sensors made by PHASEIN can be used. Oxygen exhausted oxygen sensors should be disposed according to local battery disposal stipulations.
- Never try to open the oxygen sensor device. The oxygen sensor in the IRMA probe is a disposable product, containing corrosive electrolytes and lead.
- The IRMA probe is designed only as an auxiliary means for patient evaluation. It must be used together with other vital sign and symptom evaluation equipment.
- Never place the IRMA airway adapter somewhere between the trachea catheter and the elbow; otherwise it may result in the adapter window being clogged by the patient's secretions and

operating errors.

- In order to prevent secretions and moisture from aggregating at the window and the oxygen sensor port, always place the IRMA probe at a vertical position and let the LED face upwards.
- Never use the IRMA airway adapter together with a quantitative spraying agent or spray; otherwise it may affect the light traveling of the airway adapter window.
- If an IRMA OR (without the automatic anesthetic gas identification function) user select a wrong anesthetic gas, it will result in false anesthetic gas readings.
- If the IRMA OR (without the automatic anesthetic gas identification function) is applied to a mixed gas containing several anesthetic gases, it will result in false anesthetic gas readings.
- Mobile and radio frequency communication equipment will affect measurement. You should make sure that the IRMA probe is used in the electromagnetic environment designated in this operating instruction manual.
- Never disinfect the IRMA probe or soak it into a liquid.
- The IRMA oxygen cell and the IRMA airway adapter are not germ free devices. Never carry out high pressure sterilization on the equipment; otherwise it will result in equipment damage.
- Even if an IRMA probe has not been used, do not install an oxygen exhausted oxygen cell on the probe.
- Do not stretch the sensor cable.
- Do not run this equipment beyond the temperature environment designated by this operating instruction manual.
- (USA): According to the federal law, this product can only be sold by doctors or based on prescriptions.

16.13 Airway Obstruction

When the anesthetic gas airway is obstructed, on the screen there will be such a prompt message as “The anesthetic gas airway is obstructed”; under such a circumstance, replace the Nomoline sampling tube.

Warning

- Do not use the ISA gas analyzer together with a quantitative spraying agent or pulverization treatment; otherwise it may result in the clogging of the germ filter.

16.14 Discharging Waste Gases

When nitrous oxide and/or an anesthetic gas is used, you should prevent these gases from polluting the operating room. Usually the gas discharging outlet should be connected to (via the gas discharging pipe connected to the sample gas outlet of the host equipment):

A discharging system (used for discharging collected gases) or the patient circuit (used for the back flowing of collected gases).

 **Warning**

- Anesthetics:** When an anesthetic which is being used or a patient who recently used an anesthetic is measured, the gas discharging hole on the module must be connected to a waste gas processing system or the patient circuit (on the anesthesia machine or the respirator), so as to prevent medical personnel from inhaling the anesthetic.

16.15 Consumables


The Nomoline sampling tube cannot be used repeatedly.

Replace the Nomoline sampling tube every two weeks or when “The sampling tube is clogged” is displayed (based on whichever comes first).

16.16 Safety Symbol Information

Reminding Symbol	Symbol Text, Color Code and Description Text Format	Description
	Warning: Supplementary text.	“Warning” refers to a dangerous situation which may cause personal injury or death. A warning symbol should be in conformity with: ISO 7010-W001.
	Operation Instructions	Please refer to operation instructions
	ISA operating instructions	Please refer to ISA operating instructions
	Catalog number	

	Serial number	
	Lot No.:	
	Effective period: [YYYY-MM-DD]	The equipment can no longer be run after the date nearby the symbol.
	Temperature limitation	
	Pressure limitation	
	Humidity limitation	
	No repeated usage	
	Waste electronic and electrical equipment (WEEE) directives	Electrical and electrical equipment which should be reclaimed in accordance with 2002/96/EC directives
	Containing Pb	
	Applicable to IRMAXL oxygen sensors	
 0413	European Union	In conformity with 93/42/EEC medical equipment directives when connected to medical equipment approved by PHASEIN AB
IPX4	Waterproof rating	Waterproof performance indicated by IP classification
IP44	Level of protection from water inflow and solid foreign objects indicated by IP classification	Provide protection for tools and short wire ends (>1mm). Provide protection for water spraying towards all directions.
Rx ONLY	Sold based on prescriptions	Warning (USA): According to the federal law, this product can only be sold by practicing medical personnel or based on prescriptions.

	Recommended storage temperature	Applicable to IRMA XL oxygen sensors
	Limited temperature rise	Below room temperature (23°C), the maximum surface temperature of the probe is: - IRMA CO ₂ , 41°C/106°F - IRMA OR, 45°C/113°F - IRMA AX+, 50°C/122°F
	Carbon dioxide	ISA equipment measures CO ₂ only
	Multi-gas (AX+ or OR+)	ISA equipment can measure multiple gases
	Sigma multi-gas technology	PHASEIN Sigma multi-gas technology is applied in the product
	Gas inlet	
	Gas outlet (discharge)	
	{Patient circuit connection	Illustration of the connection of Nomoline to the patient circuit
	Connected to ISA	Illustration of the connection of Nomoline to ISA
	Non germ free, no containing latex	This product is not a germ free device, does not contain latex.

16.17 Patents and Trademarks

(1) Patent Statement

PHASEIN AB owns the following patents for relevant products described in this operating instruction manual: SE519766; SE519779; SE523461; SE524086. Other patents are being applied.

(2) Trademark

PHASEIN IRMA™, PHASEIN ISA™, PHASEIN XTP™, Sigma Multigas Technology™, LEGIT™,

Nomoline™, IRMA EZ Integrator™, PHASEIN GasMaster™ and ISA MaintenanceMaster™ are trademarks of PHASEIN AB.

Tygothane® is a registered trademark of Saint-Gobain Performance Plastics Corporation.

16.18 Maintenance

The user should verify gas readings regularly; If finding any problem, please contact an engineer of the manufacturer for maintenance.

16.19 Cleaning the Analyzer

The “Plug in and measure” ISA sidestream gas analyzer should be cleaned regularly. Use ethanol or isopropyl alcohol with a maximum concentration of 70% and a wet rag to clean the analyzer..

In order to prevent the cleaning liquid and dust from entering into the ISA gas analyzer from the LEGI interface, the Nomoline sampling tube should be connected all the time during analyzer cleaning.

Before cleaning the IRMA probe, take off the disposable IRMA airway adapter.

Warning

- The Nomoline sampling tube is not a germ free device. In order to prevent the sampling tube from causing damages, please never carry out high pressure sterilization on any part of the sampling tube.
- Never sterilize the ISA sidestream gas analyzer and the IRMA probe or soak them into a liquid.
- The IRMA oxygen sensor and the IRMA airway adapter are not germ free devices. Never carry out high pressure sterilization on the equipment; otherwise it will result in equipment damage.

Chapter 17 IBP Monitoring

17.1 General

The usually invasive blood pressure monitoring can monitor arterial blood pressure, central venous pressure, pulmonary arterial pressure, left atrial pressure, right atrial pressure and intracranial pressure.

The monitor can be directly used for measuring blood vessel pressure (diastolic pressure, systolic pressure, average pressure)

Insert two IBP modules into the host monitoring equipment at the same time, and call out 4 channels of IBP pressure waveforms, as shown in the figure below:

Figure 17-1 4-Channel IBP Monitoring Screen

17.2 Understanding IBP Modules


Figure 17-2 IBP Module Interface

This monitor can support two IBP modules (4 channels) simultaneously; when two IBP modules are inserted, 4 channels of pressure waveforms and 4 groups of data can be displayed on the instrument screen.

17.3 Caveats for IBP Monitoring

Warning

- When you apply an accessory, make sure that the selected accessory meets medical instrument safety requirements.
- When connecting or using an accessory, you should avoid touching any metal part connected to an electric appliance.
- When the monitor is connected to high frequency electrosurgical equipment, in order to prevent electric leakage from burning the patient, you should prevent the sensor of the monitor from contacting a cable or high frequency electrosurgical equipment.
- Do not repeatedly use a disposable pressure sensor.
- Before starting monitoring, you should check if the sensor cable is normal. When you pull the sensor out from the receptacle, the screen will display such a prompt message as “The IBP sensor comes off”, and issue an alarm sound.
- If a liquid (not the solution used for filling the pressure tube and sensor) is spattered onto the equipment or accessory, especially when the liquid may enter into the sensor or monitor, please contact the maintenance department of the hospital.

Attention

- Only pressure sensors designated in this operating instruction manual can be used.

The designated sensor is provided with an electric shock preventing function (especially preventing a leaking electric current) and prevents the influence of a cardiac defibrillator; it can be used in surgical operations. When the patient is being defibrillated, the pressure wave may be temporarily disturbed; after defibrillation, the monitor will work normally, and the monitor operating mode and user settings will not be affected.

 Attention

- No matter a new sensor or a used sensor, it should be calibrated regularly in accordance with hospital regulations.
- When invasive blood pressure monitoring starts, the user should carry out zeroing on the transducer; during monitoring, the user should maintain the pressure sensor and the heart at the same level all the time; in order to prevent tube clogging, the user should keep injecting heparin saline to wash the tube and maintain the unobstructed condition of the pressure measurement path, and securely fixed the tube in order to prevent it from moving or coming off, which will affect invasive blood pressure measurement.

17.4 Steps of Monitoring

(1) Preparation for Measurement:

- (1) Insert the IBP module into the hoist monitoring equipment.
- (2) Insert the cable into the corresponding receptacle, check if the power supply for the monitor is turned on.
- (3) Prepare the pressure tube and sensor; the method is as follows: fill up the system with normal saline; make sure that there is no bubble in the tube system.
- (4) Connect the patient tube to the pressure tube; make sure that there is no air in the tube and the pressure tube or the sensor.

 Warning

- If there are bubbles in the pressure tube or the sensor, wash the system again with the filled liquid.

- (5) Place the sensor and the heart at the same level, approximately at middle axillary line.
- (6) Make sure that you have select a correct designation; refer to the next section for details.
- (7) Zero the sensor; refer to the next section for details.


Figure 17-3 IBP Monitoring

17.5 Setting IBP

17.5.1 Changing the Pressure Designation

- (a) Select [CH1] , [CH2] , [CH3] or [CH4] in the waveform area, in the [CH×: ×Waveform] menu popping up, select [Pressure Designation] → select a suitable designation in the following table.
- (b) Upon completion of setting, press the key “X|” to exit.

Designation	Description
ART	Arterial pressure
PA	Pulmonary arterial pressure
CVP	Central venous pressure
RAP	Right atrial pressure
LAP	Left atrial pressure
ICP	Intracranial pressure
P1/P2	Expansion pressure

17.5.2 Setting the Filtering Mode

(a) Select [CH1], [CH2], [CH3] or [CH4] in the waveform area, in the [CH×: ×Waveform] menu popping up, select [Filtering Mode] → [No Filtering], [Smooth Filtering] or [Normal Filtering].

(1) No filtering: Waveforms are not filtered; original waveforms are displayed.

(2) Smooth filtering: Smooth waveforms can be acquired.

(3) Normal filtering: Relatively smooth waveforms can be acquired.

(b) Upon completion of setting, press the key “X” to exit.

17.5.3 Setting Pressure Scales

Pressure scale adjustment: The IBP measurement range allows the user to adjust waveforms' positions on the screen by adjusting the “Upper Scale” or “Lower Scale”. Waveform scales are given in the IBP waveform area; the top and the bottom of each channel of IBP waveform represent the waveform's upper limit scale and lower limit scale, whose values are settable. The middle scale is determined by the upper scale and the lower scale; usually it is not optional.

(a) Select [CH1], [CH2], [CH3] or [CH4] in the waveform area, in the [CH×: ×Waveform] menu popping up, select [Upper Scale] and [Lower Scale] respectively → set an appropriate range.

Pressure designations of various channels: [ART], [PA], [CVP], [RAP], [LAP], [ICP], [P1] and [P2], which are all settable to a scale range: -

- Upper scale: The scale upper limit represents a pressure value. The selection range is the measurement range of the current pressure.
- Lower scale: The scale lower limit represents a pressure value. The selection range is the measurement range of the current pressure.
- Medium scale: The scale middle limit represents a pressure value. The range of this item is determined by the set range of the upper and lower scales.

(b) Upon completion of setting, press the key “X” to exit.

Attention

- The scale upper limit shall not lower than the lower limit.
- The scale lower limit shall not higher than the upper limit.

17.5.4 Setting the Pressure Unit

- Select [CH1], [CH2], [CH3] or [CH4] in the waveform area, in the [CH×: ×Waveform] menu popping up, select [mmHg] or [KPa];
- Upon completion of setting, press the key “X” to exit.

17.5.5 IBP Zeroing

- Channel 1 zeroing: When IBP1 and IBP2 modules are inserted into the receptacles, it is necessary to zero IBP.
- Channel 2 zeroing: When IBP1 and IBP2 modules are inserted into the receptacles, it is necessary to zero IBP.

17.5.6 Calibration

- Channel 1 pressure value: Numerical values adjustable in the range of 80~300.
- Channel 2 pressure value: Numerical values adjustable in the range of 80~300.
- Channel 1 zeroing
- Channel 2 zeroing

Attention

- The user should ensure that the sensor has been calibrated by zeroing before measurement; otherwise the instrument will not have an effective zero, which will result in inaccurate measurement data.

1) Sensor Zeroing:

Use the shuttle button to press the “IBP Zeroing” button, the system will start zeroing.

Zeroing Keys:

- ① Before starting zeroing, turn off the 3-way switch at the patient side.

- ② Before starting zeroing, the sensor must be exposed to the atmosphere.

The sensor must be placed at the same level as the heart, approximately at the middle axillary line.

- ③ Zeroing should be carried before monitoring begins, and at least once everyday (zeroing must be carried out every time the cable is pulled out).


Figure 17-4 IBP Calibration Connection Schematic

- 2) Mercury barometer calibration procedure: (Refer to Figure 17-4)

Mercury barometer calibration should be carried out when a new sensor is used, or based on the cycle designated by hospital regulations.

The purpose of calibration is to ensure that the system can provide accurate measurement results.

Before starting mercury barometer calibration, pressure zeroing must be carried out.

If you need to execute this procedure by yourself, you should have the following equipment.

- A standard sphygmomanometer
- A 3-way switch
- An approximately 25cm long pipe

 Warning

- **During patient monitoring, never, never, never, never perform the following operations:**
 - (1) Turn off the 3-way switch which is exposed to the atmosphere in order to carry out calibration by zeroing.
 - (2) Connect the pipe to the sphygmomanometer.
 - (3) Make sure that the connection to the patient has been disconnected.
 - (4) Connect a 3-way switch to a 3-way connector which is not connected to the patient tube (when the patient is being monitored). Connect a syringe to this 3-way connector, and connect the sphygmomanometer and pipe to the other end.
 - (5) Turn on the end to the sphygmomanometer
 - (6) In the pressure calibration menu, select the channel you want to calibrate, and adjust this channel's pressure value to be calibrated.
 - (7) Inflate so that the mercury column rises up to the pressure value set in the menu.
 - (8) Adjust repeatedly, until the value set in the menu and the pressure value shown by the sphygmomanometer are equal.
 - (9) Press the Calibrate button on the calibration menu; the instrument will start calibration.
 - (10) Wait for the calibration result. Based on the prompt message, determine a corresponding measure to be taken.
 - (11) After calibration is completed, remove the sphygmomanometer pipe and the attached 3-way connector.

17.6 Maintenance and Cleaning

 Warning

- **Before cleaning monitors or sensors, you must switch off and disconnect the AC power.**

17.6.1 Washing of Repeatedly Usable Pressure Sensors

After pressure monitoring is done, remove the pipe and cap on the sensor, rub the membrane of the sensor

with water. The sensor and cable can be scrubbed and soaked with soap water and the cleansers listed below, so as to make them clean.

Diluted ammonia water

Sodium hypochlorite (Bleaching powder used for cleansing)

Hydrogen Peroxide 3%

Alcohol 70%

Do not soak the connector into any liquid. When washing is done, let sensor completely dry before storing it. It is not necessary to consider it abnormal for the cable to have slight decolonization or for the surface to have a temporarily increased viscosity. If the tape residual on the sensor cable must be removed, using a double-sided viscose scavenger is very effective; using it carefully can minimize its damage to the cable.

The use of acetone, ethyl alcohol, ammonia water and chloroform or other strong solvents are not recommended, because they will damage the ethylene based cable as time passes.

Attention

- If what is used is a disposable sensor or cap, please never re-disinfect or reuse it.
- In order to protect the environment, disposable sensors or caps used must be reclaimed or disposed appropriately.

17.6.2 Sterilization

(1) Disinfection by Using a Chemical Solution

Remove obvious stains according to the procedure described above. Select a disinfectant which your company thinks effective for disinfecting equipment in the operating room; Buffer glutaraldehyde (e.g. glutaraldehyde or antiseptic) has been found very effective. Do not use a quadrivalentcation detergent, e.g. chlorinated benzene armour hydrocarbon amine. If the entire equipment needs to be disinfected, please soak the sensor (except electrical connectors) into a disinfectant at the recommended time. Make sure that the cap has been removed, then rinse all the components of the sensor (except electrical connectors) with germ free water or normal saline. The sensor must be totally dry before it is stored.

(2) Disinfection by Using a Gas

Disinfection by using a gas should be adopted in order to become totally germ free.

Remove obvious stains according to the cleaning procedure described above. When using a ethylene oxide gas disinfectant, in order to restrain the formation of Ethylene glycol, the sensor should be totally dry.

Do follow the operating method provided by the gas disinfectant manufacturer.

 Warning

- **The disinfectant temperature should not exceed 70°C(150°F), since the plastic in the pressure sensor will deform or melt at a temperature higher than this temperature.**

Chapter 18 ICG Monitoring

18.1 General

The basic principle of ICG (noninvasive cardiac output) measurement bases on chest impedance blood flow diagram (chest bioelectrical impedance technology: an indirect measurement method of TEB, in which cardiac ejection output per stroke is calculated by utilizing the change of chest blood flow impedance caused by cardiac ejection, and then cardiac output and other hemodynamic parameters are calculated..

The ICG module injects a high frequency AC signal into the chest of the human body, derives a chest impedance blood flow diagram via the change of impedance produced by the periodical change of blood flow.

18.2 Understanding the ICG Module



Figure 18-1 ICG Module Interface

18.3 Safety Information

- The ICG monitoring is applicable only to adults with a height ranging from 122cm to 229cm and a weight ranging from 30kg (67lb) to 159kg (341lb).
- The ICG monitor cannot be used on patients who have been installed with a minute respiratory volume type (MV) pacemaker and the MV sensor function can be turned on.
- During patient monitoring, the conductive gel on the sensors should not contact any electrically conductive substance.
- The ICG sensors can be used on one patient only.

18.4 ICG Parameters

18.4.1 Hemodynamic Values

The patient data screen displays latest hemodynamic parameters, which include the followings:

BP	Search blood pressure (SBP), diastolic blood pressure (DBP) and mean arterial pressure (MAP)	VI	Speed Index
CO	Ventricular output (_I and _Z)	ACI	Acceleration index
CI	心 Ventricular output (_I and _Z)	LCWI	Left ventricular working index (_I and _Z)
SI	Stroke index (_I and _Z)	PEP	Pre-ejection period
SV	Stroke output (_I and _Z)	LVET	Left ventricular blood ejection time
SVR	Systemic vascular resistance (_I and _Z)	PaOP	Pulmonary artery blocking pressure
SVRI	Systemic vascular resistance index (_I and _Z)	CVP	Ventricular venous blood pressure
hr	Heart rate	VEPT	Charged chest tissue volume
TFC	Chest fluid component	STR	Contraction time rate
TFI	Chest fluid index		
This HR value is not acquired by the ECG module, but directly acquired by the ICG module.			

18.4.2 List of Hemodynamic Parameters

Please refer to the contents of “Product Specifications” for the list of hemodynamic parameters.

18.5 Noninvasive Cardiac Output Display

The screen displays a parameter area for 1 channel of chest impedance waveform and ICG; the parameter area displays the values of primary parameter CI and 2 secondary parameters, as shown in the figure below:

Figure 18-2 ICG Monitoring

18.6 Affecting Factors

The following circumstances may affect ICG monitoring data:

- Infectious shock;
- Aortic regurgitation
- Severe hypertension (mean arterial pressure: > 130mmHg);
- The patient's height and weight exceed the designated ranges;
- Aortic balloon pumps insertion;
- Patient movement, including quivering;
- Signal interference caused by wire connections or power cords;
- All sorts of thoracotomy operations which cause the changes of normal chest blood or electric current.

18.7 Steps of Monitoring

- (1) Start the host monitoring equipment.

- (2) Insert the ICG plug-in module into the host monitoring equipment; the ICG module indicator illuminates, which means that the module has been successfully connected to the host monitoring equipment; otherwise, reinsert the ICG module.
- (3) Enter into the conventional screen of the host monitoring equipment, select [Exchange Waveform] to call out the “ICG” waveform and parameters which you want to monitor (if the screen has already displayed the “ICG” waveform and parameters, this step can be skipped).
- (4) Enter into the [ICG]parameter area, select [ICG Settings]to set [Patient Information Settings] in detail. Having finished setting, simply press the “” button to exit.
- (5) Connect the patient cable to the ICG module.
- (6) Having carried out comprehensive skin preparation on the patient, install the sensors on the patient.
- (7) Connect the connection wires to the sensors on the patient’s body based on the marked numbers.

18.7.1 Skin Preparation

The quality of ICG information displayed on the monitor depends on the quality of the electric signals received by the sensors. Comprehensive skin preparation carried out on the patient is necessary for the sensors to receive a high quality signal.

First select the skin areas for the sensors to be placed on, then do the followings by referring to the following procedure:

- (1) Shave the body hair in the selected areas;
- (2) Gently scrub the surface of the skin in the selected areas, so as to remove dead skin cells;
- (3) Thoroughly clean the skin at the selected positions in order to make sure that all the oily residuals, dead cells and abrasive material are completely removed. Residual abrasive particles will become the source of interference;
- (4) Before installing the sensors, please thoroughly dry the skin by rubbing.

18.7.2 Installing the Sensors

In order to acquire good signal quality and accurate data, it is very important to install the sensors at the appropriate positions; as shown in the figure, when installing the sensors:

- (1) Vertically place the neck sensors right underneath the earlobes at both side of the neck;

- (2) Place the upper chest sensors on the intersections of the xiphoid horizontal plane and the midaxillary lines at the both sides of the chest;
- (3) These two groups of sensors must be placed at directly opposite positions (180°).


Figure 18-3 Installing ICG Monitoring Sensors

18.8 Setting the Noninvasive Cardiac Output (ICG)

18.8.1 Patient Information Setting

- (a) Select [ICG] in the parameter area, select [ICG Settings], select [Patient Information Settings], input detailed patient information.
 - Set [Height] and [Weight]. The height and weight are necessary parameters for ICG measurement. When loading the ICG module, the system will automatically check these two parameters; if they are not inputted or the inputted ranges do not meet the requirements, the system will display corresponding prompts.
 - [Central Venous Pressure] and [Pulmonary Arterial Resistance Pressure]: They can be acquired by the IBP module. If they cannot be acquired from the IBP module, the user needs to input them manually. The ranges of input are $0 \sim 100 \text{ mmHg}$.

- Mean Heartbeat: [5], [10], [20], [30], [40], [50] or [60]. Since the values of ICG parameters are acquired by averaging the measurement results of multiple heartbeats, the user needs to input the mean heartbeat manually if it cannot be acquired from the ECG. The smaller this value is, the greater influence on the values of ICG parameters will be; The greater this value is, the smaller influence on the values of ICG parameters will be;
- Data Refreshing Rate: [5], [10], [20], [30], [40], [50] or [60]. The user can select one of the data as the time of automatic refreshing by the system; data options provided by the system will be refreshed based on this time.

- (b) Input other optional values.
- (c) Upon completion of setting, press the key “✕” to exit.

18.9 Malfunctions

Problem	Possible causes	Optionale Solution
There is a week or no ECG signal	The electrode is not properly placed onto the human body; the simulator battery is dead	Place the electrode onto the human body again; then press the “Submit Patient/System Data” button; replace the simulator battery, if the problem is still there, call technical service.
ECG or ICG display is suspended	Data reset	Circulate the power or click the “Submit Patient/System Data” button. If the problem is still there, call technical service.
The electrode lead condition is “Bad”	The electrode is not properly placed onto the human body, needs to be replaced.	Place the electrode onto the human body again; then press the “Submit Patient/System Data” button; use a different external patient cable; if the problem is still there, call technical service.

18.10 Maintenance

Clean the external surfaces of the patient cable and the patient electrode wire.

 **Warning**

- **The cross infection of blood products may cause virus infection. When sterilizing any equipment which may contact with a blood product, the user should handle it very cautiously.**

- Use a soft, nonabrasive rag, dip in any permitted cleanser, and then scrub all the external surfaces.
- Let all the component completely dry before using.
- Pull the power cord when cleaning.
- Use the power switch to perform the Power Off/Standby operation.

Usable Cleaning Plans:

- Mild soap water.
- Mild disinfectant (e.g. 10% calcium hypochlorite solution)

Unusable Cleaning Plans:

- Using a cleanser containing such an ingredient as isopropanol, acetone, iodine, phenol, chlorine or ammonia water.
- Carrying out high pressure sterilization or using an ultrasonic cleanser or a ethylene oxide disinfectant.
- Soaking the equipment into any cleanser or directly using any cleanser on the surfaces.
- Using an abrasive cleanser or rag.
- Wetting the connectors.

Chapter 19 C.O. Monitoring

• Machine Profile

The measurement of cardiac output (C.O.) is to measure cardiac output and other hemodynamic parameters in an invasive way by means of conventional thermodilution. The monitor can measure blood temperature, calculate cardiac output, and carry out hemodynamic calculation.

• Measurement Principle

The measurement of cardiac output is carried out as follows: Insert a floating tube into the pulmonary artery via the vein, then inject a certain amount of low temperature injecta via the floating tube; when the injecta mixes with the blood outputted by the heart, the blood temperature will change; the cardiac output is derived by measuring the curves of blood temperature changes before and after injecta injection based on the principle of heat balance.

In cardiac output measurement, the user can select a room temperature injecta or ice water injecta, inject the injecta by using the filling system or a separate syringe. The monitor can keep 6 measurements at most; if more than 6 measurements are made before measurement result editing, the earliest measurement will be deleted.

• Understanding the C.O. Module


Figure 19-1 C.O. Module Interface

• Setting the Cardiac Output (C.O.)

◆ Setting the C.O. Coefficient

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [C.O. Coefficient], input an appropriate coefficient.
 - C.O. coefficient: Set a calculation coefficient relevant to the floating tube and injecta volume. When the floating tube is replaced, this coefficient should be adjusted according the manufacturer's instruction.
- (b) Upon completion of setting, press the key “X” to exit.

◆ Setting the Injecta Temperature Source

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [C.O. Coefficient], input an appropriate coefficient.
 - C.O. coefficient: Set a calculation coefficient relevant to the floating tube and injecta volume. When the floating tube is replaced, this coefficient should be adjusted according the manufacturer's instruction.
- (b) Upon completion of setting, press the key “X” to exit.

◆ Setting the Injecta Temperature Source

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [Injecta Temperature Source] → [Automatic] or [Manual].
 - Automatic: The injecta temperature is acquired by the probe at real time.
 - Automatic: The injecta temperature can be set manually.
- (b) Upon completion of setting, press the key “X” to exit.

◆ Setting the Injecta Temperature

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [Injecta Temperature], input an appropriate coefficient.
 - Injecta temperature: 0~25.0; when [Injecta Temperature Source] is set to [Manual], the injecta temperature needs to be set manually.
- (b) Upon completion of setting, press the key “” to exit.

◆ Setting the Temperature Unit

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [Temperature Unit] → [°C] or [°F].
- (b) Upon completion of setting, press the key “” to exit.

Attention

- When the “Temperature Unit” setting is changed, the temperature unit of other monitoring parameters will change accordingly.

◆ Setting the Time Interval

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [Measurement Interval (s)] → select an appropriate interval time.
 - Measurement interval (s): 15~300s; each interval increment is 5s, the minimum time interval between two measurements; before the set time elapses, the next measurement cannot be carried out.
- (b) Upon completion of setting, press the key “” to exit.

◆ Setting the Injecta Volume

- (a) Select [C.O.] in the parameter area, in the [C.O. Settings] menu popping up, select [Injecta Volume (ml)] → select an appropriate injecta volume: 0.1~10ml.

(b) Upon completion of setting, press the key “” to exit.

• Cardiac Output Measurement

Warning

- When an accessory is used, the user should make sure that it does not contact any electrically conductive metal.
- According to the injection volume and calculation constants corresponding to the selection of the floating tube used, make sure that the calculation constants used for measurement match with the tube used. If the floating tube is changed, please input the calculation constants of the floating tube according to the floating tube operating instruction manual.

Caution

- The cable, tube and probe accessories recommended by the manufacturer should be used.

Attention

- The settings cannot be changed during C.O. measurement.

- (1) Plug the C.O. cable into the C.O. interface at the right side of the instrument.
- (2) Connect such components as the C.O. cable, the floating tube and the syringe, as shown in the figure below:


Figure 19-2 Cardiac Output Accessory Connection

- (3) Enter into the monitor's "Main Menu" or select "Patient Information" in the lower menu column, confirm the patient's height and weight.
- (4) In the monitor's "Main Menu" or select "Measurement Settings" in the lower menu column, enter into "C.O. Settings", refer to "Figure 2-1 C.O. Setting", select "Measurement Settings":
 - ◆ C.O. coefficient: Set a calculation coefficient relevant to the floating tube and the injecta volume. When the floating tube is replaced, this coefficient should be adjusted according to the manufacturer's instruction.
 - ◆ Injecta temperature source: When it is set to "Automatic", the injecta temperature is acquired by the probe at real time.
 - ◆ Injecta temperature: When the "Injecta Temperature Source" setting is set to "Manual", set the injecta temperature here manually.
- (5) In the monitor's "Main Menu" or select "Measurement Settings" in the lower menu column, enter into "C.O. Settings", refer to "Figure 2-1 C.O. Setting", select "C.O. Measurement", open the "C.O. Measurement" window:


Figure 19-3 C.O. Measurement

In the cardiac output window, the user can measure C.O. Contents displayed in the cardiac output measurement window include:

- ◆ Measurement prompt message
- ◆ C.O.: Cardiac output measurement result
- ◆ C.I.: Cardiac index
- ◆ BT: Blood temperature
- ◆ IT: Injecta temperature
- ◆ Measurement curve
- ◆ Monitoring time
- ◆ Functional buttons in the measurement window

- (6) When “New measurement can be carried out” is prompted in the prompt message area, the user can select the “Start” button, and immediately carry out fluid injection on the patient (in less than 4 seconds); a thermodilution curve is displayed in the C.O. measurement window at real time. After each measurement is completed, its measurement results will be displayed in the historical measurement window; the user should wait for a certain time before the user can carry out the next measurement by repeating this procedure.

Attention

- **During injection, 3-way switch’s floating tube connecting end opens, its injecta connecting end closes. After measurement is completed, the floating tube connecting end closes, the injecta connecting end opens, then the injecta is sucked into the injecta.**
- **If such a prompt message as “The injecta temperature is too high” appears, the inject temperature has exceeded 30°C; the users should replace the injecta, and then carry out measurement again.**
- **It is recommended that the user carry out fluid injection within 4 seconds after the user press the “Start” button for C.O. measurement.**

Moreover, in the C.O. measurement window, you can perform the following operations:

- ◆ **Start:** Start a C.O. measurement.
- ◆ **Stop:** If the measurement cannot be completed by the system after a long time, by selecting this button the user can stop the measurement currently being carried out.
- ◆ **Delete:** The user can select and delete a historical window.
- ◆ **Calculate:** Enter into the “Hemodynamic Calculation” window, in the “Input Value”, input the value which you want to calculate, press the “Calculate” button, then you can get the calculated value in the table below; refer to the following figure:

HEMODYNAMICS
✕

INPUT

C.O.	2.27	L/min	CVP	5	mmHg
HR	60	bpm	EDV	120	ml
PAWA	8	mmHg	HEIGHT	170.0	cm
ART Mean	93	mmHg	WEIGHT	70.0	kg
PA Mean	0	mmHg			

OUTPUT CALC

C.I.	1.2	L/min/m ²	BSA	1.80	m ²
SV	37.8	ml	SI	21.0	ml/m ²
SVR	3101	DS/cm ⁵	SVRI	5581	DS·m ² /cm ⁵
PVR	281	DS/cm ⁵	PVRI	505	DS·m ² /cm ⁵
LVSW	4.3	g·m	LVSWI	2.4	g·m/m ²
RVSW	2.55	g·m	RVSWI	1.42	g·m/m ²
EF	31.0	%			

SCOPE
SAVE
RECALL

Figure 19-4 Hemodynamic Calculation

**Attention**

- For hemodynamic calculation, you need to input parameters; without inputted parameters, corresponding calculation results will not be displayed.
- After the output value is calculated, you must save it before you can carry out “Calculation Review” on the record.

• Measuring the Blood Temperature

When C.O. measurement is not carried out, the user can measure the blood temperature by utilizing the

thermistor at the end of the floating tube in the pulmonary artery.

- Insert the C.O. measurement tube probe into the interface, then the blood temperature will appear.
- Set the alarm limit of the blood temperature.

• **Affecting Factors**

(一) Measurement Error Causing Factors:

- The patient moves during measurement
- Anxious patient
- Changes of the heart rate and rhythm
- Cardiac anomalies (e.g. valvular insufficiency)
- Sacculus inflation during measurement
- Inappropriate tube positioning
- Tube damage
- Inappropriate injection time
- Incorrect and inappropriate injection volume

(二) In order to acquire accurate C.O. measurements, the followings are recommended:

- The injecta temperature must be lower than the patient's blood temperature.
- Inject the fluid quickly and steadily.
- Inject the fluid at the end of expiration.
- Wait for 1 minute between two injections, so that the baseline of the blood temperature become stable again.

• Hemodynamic Calculation

◆ Input Parameters

Abbreviation	Unit	Chinese	English Full Name
HR	bpm	心率	heart rate
C.O.	L/min	心输出量	cardiac output
PAWP	mmHg	脉动脉楔压	pulmonary artery wedge pressure

Art Mean	mmHg	脉动平均压	artery mean pressure
PA Mean	mmHg	脉动脉平均压	pulmonary artery mean pressure
CVP	mmHg	中心静脉压	central venous pressure
EDV	ml	舒张期容积	end-diastolic volume
HEIGHT	cm	/	height
WEIGHT	kg	/	weight

◆ Input Parameters

Abbreviation	Unit	Chinese	English Full Name
C.I.	L/min/m ²	心脏指数	cardiac index
BSA	m ²	体表面积	body surface area
SV	ml	每搏射血量	stroke volume
SI	ml/m ²	每搏射血指数	stroke index
SVR	DS/cm ⁵	外周血管阻力	systemic vascular resistance
SVRI	DS.m ² /cm ⁵	外周血管阻力指数	systemic vascular resistance index
PVR	DS/cm ⁵	肺血管阻力	pulmonary vascular resistance
PVRI	DS.m ² /cm ⁵	肺血管阻力指数	pulmonary vascular resistance index
LVSF	g.m	左心室每搏做功	left ventricular stroke work
LVSFI	g.m/m ²	左心室每搏做功指数	left ventricular stroke work index
RVSF	g.m	右心室每搏做功	right ventricular stroke work

		功	
RVSWI	g.m/m ²	右心室每搏作 功指数	right ventricular stroke work index
EF	%	射血分数	ejection fraction

• Maintenance

Warning

- Before washing the monitor or the C.O. cable , the power must be turned off and the AC power supply must be cut off.

Cleaning

After C.O. monitoring operation is finished, remove the floating tube and temperature probe tube on the C.O. cable, scrub the C.O. cable with water. The C.O. cable can be scrubbed and cleaned by using soap water and one of the cleansers listed below.

Diluted ammonia water

Sodium hypochlorite (bleach powder used for cleansing)

Hydrogen Peroxide 3%

Alcohol 70%

Do not soak the C.O. cable into any liquid. After cleaning is finished, let the C.O. cable completely dry before storing it. It is not necessary to consider it abnormal for the C.O. cable to have slight decolonization or for the surface to have a temporarily increased viscosity. If the tape residual on the cable must be removed, using a double-sided viscose scavenger is very effective; using it carefully can minimize its damage to the cable.

Attention

- If what are used are a disposable floating tube, temperature sensor and syringe, please never re-disinfect or reuse them.
- In order to protect the environment, disposable floating tube, temperature probes and syringes must be reclaimed or disposed appropriately.

■ Sterilization

Before disinfecting, cleaning should be carried out according to the procedure described above.

Disinfecting operation may cause damage to the monitor or the C.O. cable to a certain extent. It is recommended that disinfecting operation be carried out only when it is thought necessary by your hospital maintenance plan.

Recommended disinfectants: 70% alcohol, 70% isopropanol, 2% glutaraldehyde solution.

Caution

- The user shall not use a gas (EtO) or formaldehyde for sterilization.

Chapter 20 IoC Monitoring

- **General**

IoC monitoring employs EEG of cortical electrical activity as well as analysis methods of frequency spectrum and power spectrum, to give a quantitative bispectral index, at the same time it gives comprehensive evaluation of patient's current state of consciousness combining with factors such as the quality analysis of EEG, so as to identify the depth of anesthesia.

System supports measurement of 4 parameters: IoC, EMG, ESR, SQI, provides a waveform formed by continuous detection of EEG by IoC module, referred to as IoC EEG waveforms, as well as a IoC dynamic trend, or IoC Trend.

- **Understanding IoC Module**


Fig. 0-9 IoC module connection

- **Safety Information**

For patients with anomalies of the nervous system, patients taking psychotropic drug and children under one year of age, a curious attitude should be taken in explanation of IoC values due to limited clinical experience,

Warning

- Indicators of consciousness (IoC) should not be treated as only parameters for adjusting the amount of narcotic drugs.
- Sensors and conductive parts should not have contact with other conductive parts, including contact with earth.
- In high frequency surgery, to reduce the risk of burns, IoC monitor shall not be placed between surgical site and the electro-surgical unit return electrode.
- When defibrillator is used in patients, IoC monitor cannot be placed between defibrillator and electrode pads.

• IoC Display

IoC waveform offers two display modes: the IoC ECG waveform and the IoC trend waveform.

■ IoC ECG waveform

Fig. 0-10 IoC ECG waveform

■ IoC trend waveform

Fig. 0-11 IoC trend waveform

parameters.	English	Chinese name	Meaning	Unit	Scope
IoC		Anesthetic depth index (Consciousness index)	reflects the awake level of consciousness	/	0-99
EMG	Electromyogram	EMG activity	reflects the electric power and high-frequency artifact	dB	0-100
ESR	EEG suppression rate	Suppression rate	ESR is considered to be the percentage of time in a depressed state of EEG in the last 60 seconds.	%	0-100
SQI	Signal quality index	Signal quality index	SQI value reflects the quality of signal, and provides the reliability of BIS, SEF, SR value in the last-minute period.	%	0-100

The monitored values of these 4 parameters can be affected when “impedance test in progress, impedance value exception, lead connection test in progress or lead off” occurs.

IoC	Clinical state -
99	In a sober state
80	In a calm state
60 40	In a general Anesthesia state
0	When in the depth anesthesia state, EEG waveform is straight, patient has no brain electrical activity

• **IoC Connection**

The IoC module of this monitor receives signals using Bluetooth technology. As soon as IoC monitor is opened signals can be received when IoC module is inserted into the monitor, as shown in the following figure:


Fig. 0-12 IoC connection diagram

Connection position of patients cable and electrode tablets:


patients cable colour	Electrode location Chinese (English)
Red	Center of the forehead (Center of the forehead)
Yellow	Right side of the forehead (Right side of the forehead)
Cyan	Cheek bone (Cheek bone)

Warning

- It is very important that electrode is placed in the right position in order to maintain the proper operation of the IoC.
- Electrodes can be placed at the left or right side of the scalp.
- IoC monitors use Bluetooth technology to send and receive signals, so the distance between IoC monitor and power supply adapters as well as monitors cannot exceed 10 metres, so as to avoid the loss of signal.
- IoC monitors and power supply adapters can not be placed above patient's head in order to avoid danger to patients.

• Electrode Impedance Test Results Window

Select the parameter "IoC" and enter "IoC sets", then select "electrode impedance test results" and enter this window, as shown in the following figure:

• Fig. 0-13 electrode impedance test results window

Measuring and impedance values are displayed in electrode impedance test results window :

Symbol	State or cause	solutions
--------	----------------	-----------

	electrode impedance test passed	do not require any action.
	electrode impedance test failed lead connection loss and excessive dry skin will lead overtop of impedance measurement.	check the connection between patient cable and electrode with skin, if necessary, please clean and dry the skin.
	In electrode impedance tests	No measure is required.
	Measured electrode impedance values (kΩ)	No measure is required.

• IoC Setting

◆ Set Gain

- select "IoC" in waveform area, select the "gain" in the pop-up "IoC settings" menu, → [25uV] , [50uV] , [120uV] , [250uV] or [470uV]
- Upon completion of setting, press the key “” to exit.

◆ IoC Display Mode Setting

- Select "IoC" in waveform area, then select the " IoC display mode settings " in the pop-up "IoC settings" menu→ [IoC EEG] or [IoC trend] ,
- Upon completion of setting, press the key “” to exit.

◆ Trend Time Setting

- Select "IoC" in waveform area, then select the “ trend time” in the pop-up "IoC settings" menu→ [9min] , [18min] , [60min] , [2h] , [4h] or [8h]

(b) Upon completion of setting, press the key “” to exit.

• Cleaning and Maintenance

Wipe the monitor gently with a soft damp cloth

You cannot use abrasive material in any case, so as not to cause damage to the device.

It is recommended that product correct should be taken by some qualified service personnel authorized by the manufacturer once a year.

If you have any questions, please contact manufacturer.

Attention

- If the monitor falls, damages or suffer excessive humidity or high temperature, please stop using immediately and submit it to qualified service personnel for check.

• Malfunctions

Problem	Causes	Solutions
IoC-View unable to open	check the power cord or battery	reconnect the power cord or replace the battery.
sensor alarm	● Sensor ● Battery ● No Bluetooth connection	● check the patient lead cable. ● reconnect the power cord or replace the battery. ● connect Bluetooth. ● check if the sensor is bad or the connection fails.

		● replace the damaged sensor.
index is deleted	Sensor	● sensor impedance is too high. ● If sensor impedance is greater than 15kΩ, IoC values, ESI and EMG will be deleted. ● check if the sensor is dry. ● check for proper cleaning of the skin.
signal quality index (SQI) is low	Impedance	● low impedance can cause SQI drops to 50%. ● check the electrode connection.
external signal interference	Diathermy, EMG activity	● Many reasons may cause outside signal interference: Diathermy, EMG activity. ● Make sure devices is away from the sensor.
IoC values is higher than that of the original		● Check anesthesia delivery systems. ● Due to differences between patients, some patients require a higher drug dose. ● Sufficient quantity of information may not be full to the increased stimulation.

Appendix I: Annex

The manufacturer recommends the use of the following accessories for the use of this monitor.

 **Warning**

- Please use accessory models designated by the manufacturer. Using accessories of other models may cause damage to this monitor.
- Disposable accessories can be used only once; their repeated using may result in performance deterioration or cross infection.

(一) Accessories in Standard Configuration:

No.	Accessory Number	Accessory Name	Quantity
1	040-000148-00	Earthing wire/model B/3M	1 pc.
2	009-000074-00	Power cored/national standard/3M	1 pc.
3	040-000010-00	Blood oxygen probe/adults/NELLCOR finger clipping type	1 Set
4	040-000009-00	NELLCOR blood oxygen extension wire	1 pc.
5	040-000245-00	ECG cable/12-lead/integrated snap-on type	1 pc.
6	040-000011-00	Electrode slice for adults (15 slices)	1bag
7	115-000012-00	Blood pressure tube/2M	1 pc.
8	040-000005-00	Blood pressure cuff/adults/27-35CM	1 pc.
9	040-000003-00	Body temperature probe/ body surface/ 10K	1 pc.
10	022-000001-00	4000MAH lithium battery (boxed)	1piece
11	None	4G SD card	1 pc.

(二) Optional Accessories:

(1) Standard 6-Parameter Accessories

No.	Accessory Number	Accessory Name	Quantity
1	040-000006-00	Stylus pen	1
2	040-000007-00	Printing paper/50MM×20M	1 roll
3	040-000248-00	Power adapter	1 pc.
4	040-000034-00	Blood oxygen probe/children/finger clipping type	1 pc.
5	040-000074-00	Blood oxygen probe/neonates/binding type	1 pc.
6	040-000004-00	Blood oxygen probe/children/NELLCOR (disposable)	1 pc.
7	040-000086-00	Blood oxygen probe/ neonates/binding type	1 pc.
8	040-000087-00	Blood oxygen probe/ adults, neonates/NELLCOR binding type	1 pc.
9	040-000223-00	Blood oxygen probe/ neonates/ NELLCOR (disposable)	1 pc.
10	040-000134-00	Blood oxygen probe/ adults / finger clipping type	1 pc.
11	040-000137-00	Blood oxygen probe/ adults / soft fingerstall	1 pc.
12	040-000198-00	Blood oxygen probe/ infants/MASIMO/ (disposable) 3-20KG	1 pc.
13	040-000200-00	Blood oxygen probe/ neonates /adults /MASIMO/(disposable)lighter than 3KG or heavier than 40KG (sole type)	1 pc.
14	040-000201-00	Blood oxygen probe/ children/MASIMO/ (disposable) 10-50KG	1 pc.
15	040-000202-00	Blood oxygen probe/ adults/MASIMO/ (disposable) heavier than 30KG	1 pc.
16	040-000203-00	Blood oxygen probe/ adults/MASIMO/ (repeatedly usable) heavier than 30KG	1 pc.
17	040-000232-00	Blood oxygen probe/ neonates/ adults/MASIMO/ (disposable)	1 pc.

18	040-000320-00	ECG cable /3-lead /integrated clipping type	1
19	040-000321-00	ECG cable /5-lead /integrated clipping type	1
20	040-000322-00	ECG cable /12-lead /integrated clipping type	1
21	040-000244-00	ECG cable /3-lead /integrated snap-on type	1
22	040-000247-00	ECG cable /5-lead /integrated snap-on type	1
23	040-000310-00	Electrode slices for adults /(30 slices)	1 bag
24	040-000181-00	Electrode slices for adults /(30 slices)	1 bag
25	040-000197-00	Electrode slices for neonates /(30 slices)	1 bag
26	040-000073-00	Electrode slices for neonates /(15 slices)	1 bag
27	040-000072-00	Blood pressure tube /neonates /2M	1 pc.
28	040-000141-00	Blood pressure cuff /neonates /SIZE:4*8	1 pc.
29	040-000179-00	Disposable common cuff for neonates /3-6CM	1 pc.
30	040-000180-00	Disposable common cuff for neonates /7-13CM	1 pc.
31	040-000067-00	PHILIPS silica gel type disposable cuff 1# (3.1CM ~5.7 CM)	1 pc.
32	040-000068-00	PHILIPS silica gel type disposable cuff 2# (4.3CM ~8.0 CM)	1 pc.
33	040-000069-00	PHILIPS silica gel type disposable cuff 3# (5.8CM ~10.9CM)	1 pc.
34	040-000070-00	PHILIPS silica gel type disposable cuff 4# (7.1CM ~13.1CM)	1 pc.
35	040-000085-00	Body temperature probe / body surface/ model 10K/ neonates	1 pc.
36	040-000121-00	Body temperature probe / intracavity/ 10K	1 pc.

(2) CO₂ Accessories:

RESPIRONICS Mainstream CO₂ Accessories:

No.	Code	Name	Quantity
1	099-000005-00	CO ₂ module REF:1015928	1 pc.
2	040-000021-00	Airway adapter used for a single patient REF:6063-00	1 pc.
3	040-000022-00	Airway adapter used for a single neonate REF:6312-00	1 pc.
4	040-000023-00	Cable fixing casing REF:8751-00	5 pc.

RESPIRONICS Sidestream CO₂ Accessories:

No.	Code	Name	Quantity
1	099-000004-00	CO ₂ module REF:1022054	1 pc.
2	040-000024-00	Airway adapter (with a dehumidification tube) (adults/children) REF:3473ADU-00	1 pc.
3	040-000025-00	Module fixing clip REF:1027730	1 pc.
4	040-000026-00	Airway adapter (with a dehumidification tube) (children/babies) REF:3473INF-00	1 pc.
5	040-000027-00	CO ₂ nasal tube kit (adults) REF:3468ADU-00	1 pc.
6	040-000028-00	CO ₂ nasal tube kit (babies) REF:3468INF-00	1 pc.
7	040-000029-00	CO ₂ /nasal tube kit (adults)	1 pc.

		REF:3470ADU-00	
8	040-000030-00	CO ₂ /nasal tube kit (children) REF:3470PED-00	1 pc.

PHASEIN Sidestream ISA™ CO₂ Accessories:

No.	Code	Name	Quantity
1	040-000017-00	NOMOLINE sampling tube/ CAT.NO.108210	5
2	099-000007-00	CO ₂ module /CAT.NO.800101	1 pc.
3	115-000064-00	CO ₂ accessory kit	1 pc.

PHASEIN Mainstream IRMA™ CO₂ Accessories:

No.	Code	Name	Quantity
1	099-000006-00	CO ₂ analyzer /CAT.NO.200101	1 pc.
2	040-000216-00	Airway adapter /adults /children /CAT.NO.106220	1 pc.
3	040-000217-00	Airway adapter /infants /CAT.NO.106260	1 pc.
4	115-000063-00	CO ₂ accessory kit	1 pc.

(3) PHASEIN AG Sidestream ISA™ Accessories:

No.	Code	Name	Quantity
1	040-000017-00	NOMOLINE 采样管/ CAT.NO.108210	5
2	099-000003-00	ISA™ AX+ anesthesia module/ CAT.NO.800601	1 pc.
3	115-000059-00	Anesthesia module kit	1 pc.
4	115-000065-00	ISA™ AX+ anesthesia module accessory	1 Set
5	040-000193-00	Anesthesia module extension wire accessory	1

(4) Invasive Blood Pressure (IBP) Sensor Accessories:

No.	Code	Name	Quantity
1	040-000012-00	Invasive blood pressure cable	2
2	040-000013-00	Invasive blood pressure sensor	2

(5) Noninvasive Cardiac Output (ICG) Accessories:

No.	Code	Name	Quantity
1	040-000014-00	Noninvasive cardiac output electrode slice REF:BZ-1550-50	1 bag
2	040-000015-00	Main patient cable, part NO.5551	1
3	040-000016-00	Patient lead wire, part NO.5562	2

(6) Cardiac Output (C.O.) Accessories:

No.	Code	Name	Quantity
1	040-000237-00	Cardiac output switching cable/ SPACELABS 306655-001	1
2	095-000025-00	Disposable syringe	1 Set
3	040-000238-00	Floating tube/EDWARDS/131HF7	1
4	040-000239-00	Injecta temperature probe BD/SP4042	1 pc.
5	040-000241-00	DB injecta temperature probe seat/SP5045	1 pc.

(7) IoC Anesthesia Depth Accessories:

No.	Code	Name	Quantity
1	051-000412-00	IoC module /IoC-Wiew	1 pc.
2	051-000329-00	IoC patient cable	1
3	051-000330-00	IoC dedicated electrode slice	1 bag
4	051-000331-00	IoC power adapter	1 pc.

Appendix II Product Specification

I. Monitor Type

Name	type
Classification by protection against electric shock	It belongs to the anti-defibrillator device with internal power in the I Class, the application part of the CF type and the BF type. The application part of the CF type contains ECG, invasive blood pressure, invasive cardiac output, non- invasive cardiac output and the depth detection of anesthesia, while the remaining detection parts belong to the application part of the BF type.
Classification by medical equipment management	This monitor is the III class.
Safety standard	IEC 60601-1 and GB9706.1; IEC 60601-2-27 and GB 9706.25; IEC 60601-2-4 and GB 9706.8; IEC 60601-2-26 and GB 9706.26; IEC 60601-1-2 and YY 0505; IEC 60601-2-30 and YY 0667; IEC 60601-2-49 and YY 0668; IEC 60601-1-8 and YY 0709; ISO 21647 and YY 0601; YY 0670; YY 1079.
The degree of ingress protection	Conventional equipment (the package is the IPX0 without ingress protection)
the degree of safety in the condition of flammable anesthetic gas mixed with the mixture like air, the oxygen or nitrous oxide mixture (It's NA)	The device can not be used in the case of flammable anesthetic gas mixed with air and the mixture of oxygen or nitrous oxide.
Disinfection/sterilization	Please refer to Chapter 8 to Chapter 20 for detailed information.

methods	
Operation Mode	Run the device continuously

II. Monitor Specifications

(1) Size and Weight

Name	Specifications
Size and weight	Size: 400mm（in length）×180mm（in width）× 380mm（in height）
	Weight of the complete appliance: 9 kg(battery and module are excluded)

(2) Environmental Demands

Name	Specifications	
Work environment	range of temperature	5℃～40℃
	range of Rh	< 80%
	range of atmospheric pressure	700hPa～1060hPa
Demands on the power supply voltage	power voltage	100V-240V～
	Power frequency	50Hz/60Hz
	input current	1.9A-0.8A
	direct voltage	d.c. 16.8V
	safety fuse	FUSE T 2A.

transportation	Prevent severe shock, vibration and splash from rain and snow in the process of transportation
storage	Monitors packed should be stored around -10 °C to 40 °C, no more than 80 % in the relative humidity, non-corrosive gases and ventilated room.

(3) Displayer Specification

Name	Specifications
display	17.1 inch screen, TFT display screen, color LCD
display information	8 channels waveform at most
resolution	1280×1024 pixel

(4) Record (Optional part for purchasing)

Name	Specifications
width of chart paper	50mm
valid width of recorder	48 mm
Paper speed	25/50 mm/s
settings of realtime record	3 sec, 5 sec, 8 sec, connection
tracing waveform	2 channels
alarm trigger record	yes

(5) Battery

Name	Specifications
battery Specifications	4000mAh 16.8V lithium ion battery

Duration of charging	Require at least 6 hours in the normal use.
endurance time	It can last at least 2 hours in the fully charged condition and the normal use. Moreover it can run five minutes after the first alarm triggered due to low battery power.

(6) Date Storage

Name	Specifications	
Trend data	Short-term trend	1 hour, 1 second in resolution
	Secular trend	150 hours, 1 minute in resolution
Trend chart	1200 hours	
Parameter alarm event	There are 500 parameter alarm events and parameter waveforms connected to alarm time, whose length is optional with 8s, 16s and 32s.	
analysis results of 12 leads ECG	40 sets	
NIBP measurement data	5000 sets	
Arrhythmia review	The 100 ARR events and theirs connected waveforms, whose length is optional with 8s, 16s and 32s.	
Waveform hology review	24 hours at most. The specified time is determined by the stored waveforms and theirs quantity.	

(7) ECG Specifications

Name	Specifications
------	----------------

Lead mode	12 leads (R、L、F、N、C1、C2、C3、C4、C5、C6 or RA、LA、LL、RL、V1、V2、V3、V4、V5、V6)	
Lead style	I、II、III、avR、avL、avF、V1、V2、V3、V4、V5、V6	
waveform	2 channels	
Lead mode	5 leads (R、L、F、N、C or RA、LA、LL、RL、V)	
Lead style	I、II、III、avR、avL、avF、V、	
waveform	2 channels	
Lead mode	3 leads (R、L、F or RA、LA、LL)	
Lead style	I、II、III、	
waveform	1 channel	
Overload protection	load non-damaging within 10 seconds in the IV, the power frequency and the difference module alternating voltage (p-v)	
Measurement of breath and lead and active silencing of noise	<0.1μA	
QRS wave amplitude and the range of interphase	the range of amplitude (p-v RTI)	0.5mV~5mV
	The range of width (monitor for adult)	70ms~120ms
	The range of width (monitor for neonate/ child)	40ms~120ms
	No response to the following signals	a)signal where amplitude (neonate/ child operation mode is excluded) (p-v RTI) is less than or equal to 0.15mV

		b) signal where amplitude (neonate/child operation mode is excluded) for 1mV and width for 10ms.
power frequency' voltage tolerance	$>100\mu\text{V}$ (p-v)	
Drifting tolerance	range of pyramidal wave	4mV
	range of QRS wave	0.5 mV
	width of QRS wave	100ms
	QRS wave's repetitive rate	80bpm
range of heart rate measurement and error	Range for adult	15~300bpm
	Range for neonate/ child	15~350bpm
	error	$<$ the larger one, either $\pm 10\%$ or the 5 bpm
	The widest range of measurement in the adult mode	300 bpm
	The widest range of measurement in the neonate/ child mode	350 bpm
Range of alarm limit	adult	upper limit (lower limit + 2) ~300bpm lower limit 15bpm~(upper limit-2) bpm
	neonate/ child	upper limit (lower limit + 2) ~350bpm lower limit 15bpm~(upper limit-2) bpm

Alarm resolution	no lower than the larger one of the nominal setting, either ±10% or 5 bpm		
error of the alarm limit	the maximal error should be te larger one of the nominal setting, either ±10% or the 5 bpm		
Starting time for cardiac arrest, the high heart rate alarm and the lower heart rate alarm	<10s		
Dynamic range of inputting	amplitude of inputting signal	±5mV;	
	Speed rate (RTI)	320mV/s;	
	offset voltage of direct current	-300~+300mV;	
	Changes in outputting signal	±10%;	
	display of the non-operation condition (the degree of attenuating is presented)	No declining below 50%	
input impedance	no more than 20% signal attenuating (0.67Hz~40Hz)		
system noise (p-v RTI)	<30μV		
multichannel crosstalk	<5%		
gain control and stability	Gain selection	Whole display	×2.5mm/mV, × 5.0mm/mV, × 10mm/mV, × 20mm/mV, automatic

		persistent display	10mm/mV
	gain control	support manually replacement	
	gain variation in every minute	Not more than 0.66% per min.	
	the general gain variation in an hour	Not more than $\pm 10\%$	
Selection of timing reference and accuracy	Selection of timing reference	Permanent display	25mm/s
		non-persistent display	12.5 mm/s, 25 mm/s, 50 mm/s
	The maximum error of timing reference	$\pm 10\%$.	
output display	Width of channels	30mm	
	Aspect rate	0.4s/mV	
Input signals reconstructing accuracy	The general error of system	the larger one, either $\pm 20\%$ or $\pm 100\mu V$	
	Frequency response	Sinusoidal input	0.67~40Hz (attenuating-3dB)
		Response to the wide pyramidal wave inputting 20ms	Amplitude of the wave crest attenuates from 0 to 25 Hz

	Response to the 0.3 mV's shock in its range	deviance (RTI)	Not more than 0.1mV
		slope rate (RTI)	Not more than 0.30mV/s
	Electric pole's weight factor	Not lower than ±5%	
	hysteresis effect deviating 15 mm	Not more than 0.5mm	
calibration voltage's 1mV error	±10%		
Common-mode rejection	< 1mV(p-v RTI).		
Baseline control and stability	recovery time after resetting	3s	
	Drift rate within 10s	10μV/s	
	Baseline drift with 1h	Not more than 500μV	
	Baseline drift below work temperature	Not more than 50μV/°C	
Pacemaker's pulse suppression without overshoot	Extent: from ± 2mV to ± 700mV, width: from 0.1ms to 2.0ms, the overshoot: less than 0.05α, p- the settling time: less than 5μs and the start of pulse, the rise and fall time, all of them don't exceed 100μs ; the start time of pulse goes ahead of the QRS wave's start time in or less than 40ms, that is to say, havingthe same pulse earlier than that pacing pulses in 150ms to 250ms.		
Pacemaker's pulse suppression with overshoot	the range of it is between ±2mV and±700mV and the width between 0.1ms and 2.0ms (Mthod A) cardi tach o meter's accuracy and response to arrhythmia		
Inhibition of the pacing pulse detector	Minimum input slew rate: 320mV/s RTI		

for fast ECG signal		
display capabilities of pacing pulses	Extent: from $\pm 2\text{MV}$ to $\pm 700\text{MV}$, width: from 0.5ms to 2.0ms , the maximum rise time: $100\mu\text{S}$ and the display of ECG when 100 pacing pulses appear per minute.	Not lower than 0.2mV
ST-segment detection	detecting range	$-2.0\text{mV} \text{ -- } +2.0\text{mV}$
	detecting accuracy	in the range of $-0.8\text{ mV} \sim +0.8\text{mV}$, the detecting error is the larger one, either $\pm 0.02\text{mV}$ or $\pm 10\%$
resolution	0.01 mV	
arrhythmia types	asystole, ventricular fibrillation / tachycardia, a single ventricular ectopic, two ventricular ectopic, multiple ventricular ectopic, ventricular ectopic bigeminal heartrate, ventricular ectopic triplet heartrate, the R- on -T (R ON T) , missed beat, tachycardia, bradycardia, pacemaker not being captured and not being paced	
leakage current	$< 10\text{ uA}$	
electrosurgery interference suppression	The change of heart rate is not more than $\pm 10\%$ before the interference	
Electrotome protection	incision mode: 300W ; solidity mode: 100w , restoring time: $\leq 10\text{s}$	

heart rate algorithm	
Suppression capacity of high and big T	When tested according to section 4.1.2.1 c of the YY 1079, the cardiometer will suppress all the QRS wave group whose volatility is less than

wave	1.0mV and T wave which has 180ms interval and 350ms QT interval.
average heart rate	According to the requirements of 4.1.2.1 d part of the YY 1079, calculate the average heart rate by the following methods: if the last consecutive three RR have intervals greater than 1200 ms, get the mean from four of the most recent RR intervals to calculate the heart rate; Otherwise, it would take 12 the most recent RR intervals to cut the maximum and minimum, and get the mean to calculate the heart rate.
Accuracy of the cardiograph and response to arrhythmia	Required by 4.1.2.1 e- part of the YY 1079, after the stabilization period of 20 seconds, the heart rate values are: Figure 3 a) (bigeminal rhythm) $80 \pm 1\text{bpm}$; Figure 3 b) (slowly varying bigeminal rhythm) $60 \pm 1\text{bpm}$; Figure 3 c) (the fast bigeminal rhythm) $120 \pm 1\text{bpm}$; Figure 3 d) (two-way contraction) $90 \pm 2\text{bpm}$.
Time for response to heart rate changes	Required by 4.1.2.1 f- part of the YY 1079, not more than 10s it rises from 80 to 120 bpm; it doesn't need 10s to decline from 80 to 40bpm.
starting time of the tachycardia alarm	Required by 4.1.2.1 g part of the YY 1079, the waveform: Figure 4 a) 1 - scope: 10s; Figure 4 a) 0.5- scope: 10s; Figure 4 a) 2 - scope: 10s; Figure 4 b) 1 - scope: 10s; Figure 4 b) 0.5 - scope: 10s; Figure 4 b) 2 - scope: 10s.

(8) Breathing Specifications

Name	Specifications
------	----------------

Name	Specifications		
manner	way of Thoracic impedance		
Detection and accuracy of breathing rate	Range of detection	adult	6bpm-120bpm
		neonate/child	6bpm-150bpm
	Accuracy of detection	± 1 bpm	
accuracy and error of the alarm preset in breathing rate	adult	upper limit not narrower than	(lower limit + 2) ~ 120bpm
		lower limit not narrower than	6bpm ~ (upper limit - 2) bpm
	neonate/child	upper limit is not narrower than	(lower limit + 2) ~ 150bpm
		lower limit is not narrower than	6bpm ~ (upper limit - 2) bpm
	error	± 1 bpm	
range and error of the suffocation alarm	range	10s~40s	
	error	± 5 s.	

(9) NIBP Specifications

Name	Specifications		
Way of measurement	self-oscillation method		
range and accuracy of measurement	range of measurement for adult	systolic pressure	5.3-36kPa (40-270mmHg)
		diastolic pressure	1.3-28.7kPa (10-215mmHg)
		mean pressure	2.7-31.3kPa (20-235mmHg)

	range of measurement for child	systolic pressure	5.3-26.7kPa (40-200mmHg)
		diastolic pressure	1.3-20kPa (10-150mmHg)
		mean pressure	2.7-22kPa (20-165mmHg)
	range of measurement for neonate	systolic pressure	5.3-20kPa (40-150mmHg)
		diastolic pressure	1.3-13.3kPa (10-100mmHg)
		mean pressure	2.7-14.7kPa (20-110mmHg)
	accuracy of measurement	±5mmHg, when the non-invasive blood pressure is beyond the range, the monitor still displays properly, but does not consider accuracy.	
Range and permissible deviation of overvoltage protection	Adult mode	300mmHg	
	child mode	240mmHg	
	Neonate mode	150mmHg	
	Permissible deviation	±3mmHg	
Range and error of the alarm preset	adult	systolic pressure	upper limit (lower limit + 2) – 270bpm lower limit 40 – (upper limit – 2)mmHg
		diastolic pressure	upper limit (lower limit + 2) – 215bpm lower limit 10 – (upper limit –

			2)mmHg
		mean pressure	upper limit (lower limit + 2) – 235bpm lower limit 20 – (upper limit – 2)mmHg
	Child	systolic pressure	upper limit (lower limit + 2) – 200bpm lower limit 40 – (upper limit – 2)mmHg
		diastolic pressure	upper limit (lower limit + 2) – 150bpm lower limit 10 – (upper limit – 2)mmHg
		mean pressure	upper limit (lower limit + 2) – 165bpm lower limit 20 – (upper limit – 2)mmHg
	Newborn	systolic pressure	upper limit (lower limit + 2) – 135bpm lower limit 40 – (upper limit – 2)mmHg
		diastolic pressure	upper limit (lower limit + 2) – 100bpm lower limit 10 – (upper limit – 2)mmHg
		mean pressure	upper limit (lower limit + 2) – 110bpm

		lower limit 20 — (upper limit — 2)mmHg
	error	the larger one, either $\pm 0.1\text{kPa}$ or $\pm 1\text{mmHg}$
Mode of blood pressure measurement	Manual and automatic (periodic), continuous	
	test interval of the manual mode	1、 2、 3、 4、 5、 10、 15、 30、 60、 90、 120、 180、 240、 480min
	Continuous mode	5min

(10) SpO₂ Specifications

Name	Specifications	
range of display	1%~100%	
resolution of display	1%	
Accuracy measurement	In the scope of 90%~100%, the measurement error is $\pm 1\%$ (no motion state). In the scope of 70%~90%, the measurement error is $\pm 1\%$ (no motion state).	
Limit and accuracy of alarm preset	Upper limit of alarm	(lower limit+1)%~100%
	Lower limit of alarm	0%~(lower limit-1)%
	accuracy	$\pm 1\%$

(11) Pulse Rate Speculations

Name	Specifications
------	----------------

range and accuracy of measurement	Range	It shall be 20bpm~300bpm.
	accuracy	±1bpm
Alarm se and accuracy of pulse rate	alarm upper limit: (lower limit + 1)~300bmp	
	alarm lower limit: 20~(lower limit - 1)bpm	
	±2 bpm	

(12) TEMP Specifications

Name	Specifications	
range and accuracy of measurement for	range of measurement	0℃~50℃
	error of measurement	±0.1℃
set and accuracy of alarm	Range of alarm set	T1/T2 alarm upper limit: (lowerlimit + 1)~50.0℃ T1/T2 alarm lower limit: 0℃~(lower limit - 1)0} TD height limit: 0~50 ℃
	Alarm error	±0.1℃
Display resolution	0.1℃	
Number of channel	Double channels	

(13) CO₂ Specifications

Name	Specifications
CO ₂ 's range of measurement	0mmHg~150mmHg, 0%~19.7%, 0kPa~20kPa (at760mmHg). And air pressure is supplied by the mainframe.

CO ₂ 's resolution	1mmHg or 0.1kPa or 0.1%
CO ₂ 's accuracy	± 2 mmHg between 0mmHg ~40mmHg ; $\pm 5\%$ between 41mmHg ~70mmHg; $\pm 8\%$ between 71mmHg ~100mmHg; $\pm 10\%$ between 101mmHg~150mmHg
AwRR's range of measurement	0~150 rpm
AwRR's accuracy of measurement	1rpm
Range of alarm set	EtCO ₂ 's upper limit: (lower limit +2)~150mmHg lower limit: 0~(upper limit -2)mmHg awRR 's upper limit: (lower limit +2)~150rpm lower limit: 0~(upper limit -2)rpm INS's upper limit: 0~99 mmHg
Error of alarm set	± 0.1 kPa or ± 1 mmHg

(14) IBP Specifications

Name	Specifications	
Number of IBP channel	4 channels	
Names of channel pressure	Artery pressure, pulmonary artery pressure, central vein pressure, right atrial pressure, left atrial pressure, intracranial pressure, plus pressure.	
Range and accuracy of measurement	ART	0~40kPa (0~300mmHg)
	PA	-0.8~16kPa(-6~120 mmHg)
	CVP	-1.3~5.3kPa(-10~40mmHg)

	RAP	-1.3~5.3kPa(-10~40mmHg)
	LAP	-1.3~5.3kPa(-10~40mmHg)
	ICP	-1.3~5.3kPa (-10~40mmHg)
	P1, P2	-6.6~40kPa(-50~300mmHg)
Range and accuracy of IBP pressure measurement	±1mmHg or ±2%, take the larger one (sensor error excluded)	
Range of IPB's alarm set	ART	upper limit (lower limit +2)~300 mmHg lower limit 0~(upper limit -2)mmHg
	PA	upper limit (lower limit +2)~120mmHg lower limit -6~(upper limit -2)mmHg
	CVP	upper limit (lower limit +2)~40mmHg lower limit -10~(upper limit -2)mmHg
	RAP	upper limit (lower limit +2)~40mmHg lower limit -10~(upper limit -2)mmHg
	LAP	upper limit (lower limit +2)~40mmHg lower limit -10~(upper limit -2)mmHg
	ICP	upper limit (lower limit +2)~40mmHg lower limit -10~(upper limit -2)mmHg
	P1	upper limit (lower limit +2)~300mmHg lower limit -50~(upper limit -2)mmHg
	P2	upper limit (lower limit +2)~300mmHg lower limit -50~(upper limit -2)mmHg
error of IPB's alarm set	±0.1 or ±1mmHg of the set value	

pressure sensor	sensitivity: 5μV/V/mmHg
	range of impedance: 300~3000Ω

(15) AG Specifications

Name	Specifications		
AG measurement	feature of infrared radiation absorption		
AG preheat time	< 20 sec		
Range and accuracy of AG measurement	the following standard is suitable for a dry gas in 22 ± 5 °C and 1013 ± 40 hPa:		
	CO2	0 to 114mmHg	± (1.52 mmHg + 2% of the reading)
		114 to 190 mmHg	unspecified
	N2O	0 to 100 %	± (2 % + 2% of the reading)
	Halothane enflurane isoflurane	0 to 8 % 8 to 25 %	± (0.15 % + 5% of the reading) unspecified
	sevoflurane	0 to 10 % 10 to 25 %	± (0.15 % + + 5% of the reading) unspecified
	desflurane	0 to 22 % 22 to 25 %	± (0.15 % + + 5% of the reading) unspecified
	O2	0 to 100 %	± (1 % + + 2% of the reading)
	Breathing rate	0 to 254rpm	±1 rpm
AG resolution	CO2: 1mmHg awRR: 1rpm		
AG gas	CO2 , O2, N2O and one of five anesthetics (Enflurane, isoflurane,		

	sevoflurane, halothane, desflurane)			
range and accuracy of AG alarm set	Alarm specification	range	step	accuracy
	EtCo2	upper limit (lower limit + 2)~76mmHg lower limit 0~ (upper limit - 2)mmHg	1mmHg	±1mmHg
	FiCo2	upper limit (lower limit + 2)~76mmHg lower limit 0~ (upper limit - 2)mmHg	1mmHg	±1mmHg
	AwRR	upper limit (lower limit + 2)~100rpm lower limit 0~ (upper limit - 2)rpm	1rpm	±1rpm
	EtO2	upper limit (lower limit + 2) 20~100% lower limit 18~ (upper limit - 2)%	1%	±1%
	FiO2	upper limit (lower limit + 2)~100% lower limit 18~ (upper limit - 2)%	1%	±1%
	EtN2O	upper limit (lower limit + 2) 20~100% lower limit 0~ (upper limit - 2)%	1%	±1%
	FiN2O	upper limit (lower limit + 2)	1%	±1%

		20~100%		
		lower limit 0~ (upper limit — 2)%		
	EtHal/EtEnf/EtI So/EtSev/EtDes	upper limit (lower limit + 0.2) ~25% lower limit 0~ (upper limit — 0.2)%	0.1%	±0.1%
	FiHal/FiEnf/FiI So/FiSev/FiDes v	upper limit (lower limit + 0.2) ~25% lower limit 0~ (upper limit — 0.2)%	0.1%	±0.1%

(16) ICG Specifications

Name	Specifications		
ICG way of measurement	indirect measurement by the impedance cardiogram		
ICG range of measurement	SV: 5~250 ml/ beat HR: 40~250 bpm C.O.: 1.4~15L/min		
ICG accuracy	SV: no definition HR: ±2bpm C.O: no definition		
ICG range of alarm set	Alarm specification	range	step
	CI	upper limit (lower limit +0.1) ~ 15.0 L/min/m2 lower limit 0~ (upper limit —0.1)	01. L/min/m2

		L/min/m2	
	TFC	upper limit (lower limit + 0.1) ~150/KΩ lower limit 10~ (upper limit - 1) KΩ	1 KΩ
ICG accuracy of alarm	CI	±0.1 L/min/m2	
	TFC	±1 KΩ	

(17) C.O. Specifications

Name	Specifications	
C.O. range of measurement	C.O.:	0.1~20L/min
	BT:	25~43℃
	IT:	0~25℃
resolution	C.O.:	0.1L/min
	BT、IT:	0.1℃
accuracy	C.O.:	take the larger one, either±5% or±0.1 L/min
	BT、IT:	±0.1℃ (sensor are excluded)
C.O. the upper and lower limit of alarm	BT upper limit	(lower limit +0.1) ~43℃
	BT lower limit	23.0~(upper limit -0.1) ℃
	Step	0.1℃
accuracy of alarm	±0.1℃	

(18) IOC Specifications

Name	Specifications		
IOC range of measurement	IOC: 0-99 SQ I:0-100% EMG: 0-100dB ESR: 0-100%		
IOC the upper and lower limit of alarm	Alarm specification	range	step
	IOC	upper limit : (lower limit +0.1) ~99 lower limit: 0~(upper limit -2)	1

(19) Wireless Network

Standards	IEEE802.11g, IEEE802.11b, IEEE 802.3 compatible with Wi-Fi and WPA
range of frequency	2.4GHz~2.4835GHz
type of modulation	OFDM、BPSK、QPSK、16 QAM、64 QAM 、CCK、DQPSK、DBPSK
signal path	1-11 (America, Canada), 1-13 (China, Europe), 1-14 (Japan)

Appendix III System Alarm Information

Some of the most important physiological and technical alarm information is listed in this section, however, some might not be included.

XX represents some module name or physiological parameters in the system such as HR,ST1,ST2,PVCs,RR,TEMP(including TEMP1,TEMP2,TD),SpO2,PR,CO2(including AwRR,INS),NIBP(NS,NM,ND),AG(including CO2,O2,N2O,AA,HAL,ENF,ISO,SEV,DES),C.O.,IBP,ICG,Io,and so on.

A corresponding remedy is offered following each warning. If the problem still exists after the remedy, please contact the maintenance crew.

(一)Physiological alarm

Origin	Alarm	Alarm Levels	Causes	Solution
XX	Too high XX	An user's option	XX is beyond the upper value of the alarm limit.	Check the alarm limiting values, or the patient's current condition.
	Too low XX	An user's option	XX is below the lower value of the alarm limit.	
ECG	Too weak ECG signals	High	The patient's ECG signals are too weak for ECG analysis by the system.	Check the connections of the electrodes and the lead cords, and the patient's current condition.
	ASYSTOLE (cardiac arrest)	High	Asystole arrhythmias occurs in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	VFIB/VTAC (ventricular fibrillation/	High	Ventricular fibrillation or tachycardia occurs in the patient.	Check the patient's current condition and the connections of the electrodes and the lead

Origin	Alarm	Alarm Levels	Causes	Solution
	tachycardia)			cords.
	COUPLET (couplet premature ventricular contractions)	Average	Couplet premature ventricular contractions occur in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	BIGEMINY (bigeminy coupled rhythm)	Average	Premature ventricular contractions in a bigeminy coupled rhythm occur in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	TRIGEMIN Y (trigeminy rhythm)	Average	Premature ventricular contractions in a trigeminy rhythm occur in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	R ON T	Average	R ON T arrhythmia occurs in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	PVC (premature ventricular contractions)	Average	Premature ventricular contractions occur in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	TACHY (tachycardia)	Average	Tachycardia occurs in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	BRADY (bradycardia)	Average	Bradycardia occurs in the patient.	Check the patient's current condition and the connections of the electrodes and the lead

Origin	Alarm	Alarm Levels	Causes	Solution
				cords.
	VT > 2 (multiple premature ventricular contractions)	Average	Multiple premature ventricular contractions occur in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	MISSED BEATS (missed beats)	Average	Missed beats occur in the patient.	Check the patient's current condition and the connections of the electrodes and the lead cords.
	PNP (Pacemaker non-pacing)	Average	The pacemaker does not pace.	Check the connections of the pacemaker, the electrodes and the lead cords as well as the patient's current condition.
	PNC (Pacemaker non-capturing)	Average	The pacemaker signals are not captured.	Check the patient's current condition and the connections of the electrodes and the lead cords.
SpO2	Pulse missing	High	The patient's pulse signals are too weak for pulse analysis by the system.	Check the connection of the sensor and the patient's current condition.
RESP	RESP asphyxia	High	The patient's respiration signals are too weak for RESP signal analysis by the system.	Check the connections of the lead cords and the patient's current condition.
CO2	RESP asphyxia	High	The patient's respiration signals are too weak for RESP signal analysis by	Check the connections of the lead cords and the patient's current condition.

Origin	Alarm	Alarm Levels	Causes	Solution
			the system.	
Others:	Alarm pause	None	Alarming is artificially suspended.	Press the Silence/Reset button again to restore the normal alarming function.

(二) Technical Alarm Information

Origin	Alarm	Alarm Levels	Causes	Solution
XX	Initialization errors of the XX module	High	X error occurs in XX module initialization.	Restart the machine for another try. If the same problem still exists, please contact the manufacturer.
	XX module communication failure	High	The XX module cannot communicate with the main system.	
	XX communication errors	High	The XX module cannot communicate normally with the main system.	
XX	XX alarm limit errors	High	The parameters of the XX alarm limits are altered accidentally.	Please contact the manufacturer for maintenance.
XX	XX out-of-range measurement	High	The XX parameter goes beyond the measurement range of the system.	Please contact the manufacturer for maintenance.
ECG	ECG lead shedding	Low	The ECG lead cords do not have a good	Check the connections of the ECG lead cords.

Origin	Alarm	Alarm Levels	Causes	Solution
			connection.	
	ECG YY lead shedding (YY represents V, LL, LA, RA, I, V2, V3, V4, V5 and V6)	Low	The ECG YY lead cords do not have a good connection.	Check the connections of the YY lead cords.
	Too strong interference in ECG	High	Strong interference appears in the ECG signals.	Check ECG lead wire connections and the current status of patients, whether there is a large action.
	ARR on-going learning	None	The formation of the ORS wave template is in process.	The status will terminates after the ARR learning finishes.
	ECG lead shedding	Low	The ECG lead cords do not have a good connection.	Check the connections of the ECG lead cords.
SpO2	SpO2 shedding off the finger	Low	The SpO2 sensor has fallen off the finger.	Check the connection of the SpO2 sensor.
	No probe connected to SpO2	Low	SpO2 sensor is not connected well.	Check the connection of the SpO2 sensor.

Origin	Alarm	Alarm Levels	Causes	Solution
	SpO2 sensor shedding	Low		
	Reset due to NELL breakdown	Low	The system is reset due to NELLCOR module breakdowns.	If the system cannot be reset or the same problem still exists after system re-start, please contact the manufacturer for maintenance.
	Pulse searching	None	The SpO2 sensor does not have a good connection or the patient has moved their arm.	Check the connection of the SpO2 sensor and the patient's current condition.
TEMP2	TEMP1 sensor shedding	Low Low	The TEMP sensor does not have a good connection.	Check the connection of the TEMP sensor.
	TEMP2 sensor shedding	=		
NIBP	NIBP self-checking	Low	NIBP errors occur during initialization.	Select the Reset function in the NIBP menu. If the same problem still exists after reset, please contact the manufacturer for maintenance.
	NIBP communication errors	Low	NIBP communication errors occur.	Select the Reset function in the NIBP menu. If the same problem still exists after reset, please contact the

Origin	Alarm	Alarm Levels	Causes	Solution
				manufacturer for maintenance.
	Cuff loosening or connection failure	Low	The NIBP cuff does not have a good connection.	Please reconnect the NIBP cuff.
	Gas leak	Low	The NIBP cuff has a poor connection, or there are gas leaks in the gas path.	Check the connections of different sections, or replace the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	Air pressure errors	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis calculation.	Check the connections of different sections, or replace the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	Too weak signals	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis calculation.	Check the patient type setup, the connections of different sections, or replace the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	Out-of-range pressure	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis	Check the connections of different sections, or replace the cuff. If the problem still exists, please contact the manufacturer for maintenance.

Origin	Alarm	Alarm Levels	Causes	Solution
			calculation.	
	Arm move	Low	The patient has moved their arm.	Check the connections of different sections, or place the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	Overpressure protection	Low	There may be folding in the gas path.	Check the gas path and the patient's current condition, and then perform measurement again. If the alarm still exists, please contact the manufacturer for maintenance.
	Signal saturation	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis calculation.	Check the connections of different sections, or place the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	Pump gas leak	Low	There is gas leak in the NIBP gas path.	Check the connections of different sections, or replace the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	NIBP system failures	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis	Check the connections of different sections, or place the cuff. If the problem still exists, please contact the

Origin	Alarm	Alarm Levels	Causes	Solution
			calculation.	manufacturer for maintenance.
	Measurement timeout	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis calculation.	Check the connections of different sections, or place the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	Cuff model errors	Low	This may be caused by the mismatching between the applied cuff and the patient type which has been set up in the system.	Check the patient type setup, the connections of different sections, or replace the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	NIBP measurement failures	Low	Errors occur during curve measurement, which incapacitates the system in performing measurement analysis calculation.	Check the connections of different sections, or place the cuff. If the problem still exists, please contact the manufacturer for maintenance.
	NIBP mistaken reset	Low	Unlawful reset occurs during NIBP measurement.	Check the NIBP gas path to make sure whether there is a blockage or not, and then perform measurement again. If the problem still exists, please contact the manufacturer for maintenance.
	NIBP	Low	NIBP errors occur during	Select the Reset function in the NIBP menu. If the same

Origin	Alarm	Alarm Levels	Causes	Solution
	self-checking		initialization.	problem still exists after reset, please contact the manufacturer for maintenance.
AG	starting calibration	Low	AG analyzer is in calibration.	Waiting for the end of calibration.
	indoor air calibration is needed.	Low	oxygen sensors require calibration in the air	oxygen sensors need to be calibrated in the 21% air.
	oxygen range calibration is needed.	Low	oxygen sensor need to be calibrated in pure oxygen .	oxygen sensors calibration need to be in 100% pure oxygen, if the customer does not have the conditions, it can be return to the factory for maintenance.
	Oxygen sensor error	Low	module fails	Reseat power module or restart, if the fault persists, please contact the manufacturer for maintenance.
	sampling tube blocking	Low	sampling tube blocking	check and replace the sampling tube, if the fault persists, please contact the manufacturer for maintenance.
	replace adapter	Low	adapter abnormal	check and replace the adapter, if fault persists, please contact the manufacturer for maintenance.。
	no sampling	Low	sampling tube is not connected or in bad	checks and connect the sampling tube, if the fault

Origin	Alarm	Alarm Levels	Causes	Solution
	tube		contact.	persists, contact the manufacturer for maintenance.
	adapter-free	Low	adapter is not connected or in bad contact.	checks and connect the adapter, if the fault persists, contact the manufacturer for maintenance.
	oxygen port failure	Low	module fails	Reseat power module or restart, if the fault persists, please contact the manufacturer for maintenance.
	CO2 exceed precision scope	Low	the measured value exceeds the claimed scope of accuracy	Please follow the precision range of measurement claimed by manufacturers
	O2 exceed precision scope	Low		
	N2O exceed precision scope	Low		
	AXsuper precision scope	Low		
	Inside Temperature exceed scope	Low	module fails	returns to factory for maintenance.
	atmospheric pressure	Low	module fails	returns to factory for

Origin	Alarm	Alarm Levels	Causes	Solution
	super the scope			maintenance.
	requires zero	Low	AG Analyzer requires zero	enter AG settings for zero Here zero is equal to calibration.
	anesthetic concentration identification is uncertain	Low	AG Analyzer does not recognize the existence of other gas	check the gas.
	replace oxygen sensor	Low	beyond usage time	replace oxygen sensor.
	software errors, restart the sensor	Low	software error occurs	Reboot
	hardware errors, update sensor	Low	module fails	Reseat oxygen sensor or reboot.
	motor velocity overrun	Low	module fails	return to factory for maintenance.
	factory calibration lost	Low	module fails	return to factory for maintenance.
	2 kinds of anesthetic gas are detected	Low	There is only one main anesthetic gas waveform and value shown in the interface. AG detects 2	check the patient's physiological signs and anesthesia.

Origin	Alarm	Alarm Levels	Causes	Solution
			kinds of anesthesia gas for tips.	
	Anesthesia module is on standby	Low	AG Analyzer operating modes is on standby.	measurement mode of AG Analyzer is selected.
C.O.	BT head off	Low	BT probes is not connected or in bad contact.	check and connect BT probe, if fault persists, contact the manufacturer for maintenance.
	IT probes off	Low	IT probes is without connection or in bad contact.	Check and connect IT probes, if fault persists, contact the manufacturer for maintenance.
	injection temperature exceed scope	Average	injection temperature is greater than 30 ° c	Cool the injection, for example, put a container with ice water below the injection duct to cool.
	injection temperature timeout	Average	injection took too much time since "start" too long.	injections time should not exceed 4 seconds since the "start" measurement.
	baseline drift	Average	It is likely to be interfered by high frequency equipments from electrosurgery.	monitor and C.O. sensor do not near electrosurgical equipments .
	baseline noise	Average	It is likely to be interfered by high frequency equipments from electrosurgery.	monitor and C.O. sensor do not near electrosurgical equipments .
IBP	IBP1 successfully zero	Low	IBP transducer zero success	None

Origin	Alarm	Alarm Levels	Causes	Solution
	IBP2 successfully zero			
	IBP1 zero fails	Average	IBP transducer zero failed	Recalibrate.
	IBP1 zero fails			
	IBP1 zero	Low	IBP sensor is being calibrated	Waiting for the calibration end
	IBP2 zero			
	IBP1 lead off	Low	IBP sensor is in bad contact.	checks or reseal IBP sensor, if the fault persists, contact the manufacturer for maintenance.
	IBP2 lead off			
ICG	ICG lead off	Low	ICG sensor is not connected or in bad contact.	check and connect ICG sensor, if the fault persists, contact the manufacturer for maintenance.
Other alarm information	Button errors	High	A system fault occurs.	Restart the monitor. If the problem still exists, please contact the manufacturer for maintenance.
	Too high voltage at 5 v	High	A power fault occurs in the system.	Restart the monitor. If the problem still exists, please contact the manufacturer for maintenance.
	Too low voltage at 5 v	High		

Origin	Alarm	Alarm Levels	Causes	Solution
	Too high voltage at 12 v	High		
	Too low voltage at 12 v	High		
	Too low battery energy	High	The electric energy of the batter is low.	Connect to alternating current and recharge the battery.. If the problem still exists after 6-hour recharging, please contact the manufacturer for maintenance.
	The battery can work for 20mins	Average		
	10-minute battery-powered endurance	High		
	Low battery energy	High		
	Automatic system shutdown: xxxS	High		

(三)System Prompts

Origin	Alarm	Alarm Levels	Causes
NIBP	Manual measurement...	None	System prompt information is only prompted for a certain function or an on-going operational procedure.
	Alignment	None	
	Leak detection	None	
	Manual reset...	None	
	Consecutive measurement...	None	
	Please press the Restart button	None	
	Error reset	None	
	Modular reset...	None	
	Automatic measurement...	None	
	Adjustment termination	None	
	Gas leak detection termination	None	
	Measurement termination	None	
	Overpressure protection	None	
	Reset failure	None	
CO ₂ calibration	Please wait during automatic zero adjustment x S	None	
C.O.	measuring	Low	
Alarm	Keyboard alarm disabled	None	

Origin	Alarm	Alarm Levels	Causes
disabled prompt information	ECG alarm disabled	None	
	Monitoring prohibition during adjustment	High	
	RESP alarm disabled	None	
	TEMP alarm disabled	None	
	SpO2 alarm disabled	None	
	NIBP alarm disabled	None	
	CO2 alarm disabled	None	
	Completely silent alarming	None	
Module loading prompt information	Key set module loading	None	
	ECG module loading	None	
	SpO2 module loading	None	
	NIBP module loading	None	
Module unloading prompt information	Module unloading	None	
Other prompt information	Demonstration	None	
	The screen has been locked. Press any key to unlock!	None	

Origin	Alarm	Alarm Levels	Causes
	The screen has been unlocked!	None	
	System booting time	None	
	Parameter module reset	None	

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Юридический адрес: South Side of Floor 7 Bldg. 5, Floor 6 Bldg. 4, Fourth Industrial Zone of Nanyou, Nanshan Avenue, Nanshan District, Shenzhen City

Фактический адрес: Floor 6 Bldg. 4, Fourth Industrial Zone of Nanyou, Nanshan Avenue, Nanshan District, Shenzhen City

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