STORZ

KARL STORZ — ENDOSKOPE

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**1 Important information
for users of KARL STORZ
devices**

Thank you very much for your expression of trust in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern high-quality piece of equipment from KARL STORZ.

This instruction manual is intended as an aid in setting up, connecting, and operating the UNIDRIVE® S III ENT motor system with the various application instruments

- INTRA drill handpiece (2524xx, 2525xx)*
- Stammberger-CASTELNUOVO DrillCut-X shaver handpiece (40711050) with shaver attachment (412xx, 413xx)
- Stammberger-CASTELNUOVO DrillCut-X shaver handpiece (40711050) with sinus drill attachment (413xx)
- Drill Cut-X II shaver handpiece (40712050/40712055) with shaver attachment (412xx, 413xx)
- Drill Cut-X II shaver handpiece (40712050/40712055) with sinus drill attachment (413xx)
- Stammberger-Sachse intranasal drill (660000)*
- Micro saw (254000)*
- Micro sagittal saw (254100)*
- Osseoscalpel (254200)*
- Micro compass saw (254300)*
- Dermatomy (253xxx)*

(*in conjunction with the EC Micromotor 20711033).

Please read these instructions carefully. Keep this manual in a convenient and conspicuous location and in its protective casing close to the device.

Users of the various application instruments should familiarize themselves with the operating instructions included with the respective devices.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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
**3 Controls, displays,
connectors, and their uses**

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Controls, displays, connectors, and their uses

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- ① Power switch for turning the unit on and off
- ② TFT touch screen
- ③ Connection socket for motor/handpiece M1
- ④ Connection socket for motor/handpiece M2
- ⑤ Irrigation pump
- ⑥ Holder for irrigator rod
- ⑦ Connection socket for footswitch
- ⑧ USB port (for maintenance only)
- ⑨ Potential equalization connection
- ⑩ Power socket
- ⑪ Fuse holder
- ⑫ SCB® connection sockets*
- ⑬ Dual-pedal footswitch
- ⑭ Button for changing direction of rotation

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*** NOTE:** The KARL STORZ SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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Read the instructions carefully before operating
equipment!

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Electronic information product pollution control
(China RoHS)

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DANGER: Risk of explosion if used in the
presence of flammable anesthetics.

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CAUTION: To reduce the risk of electrical
shock, do not remove cover.

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Refer service work to qualified service per-
sonnel.

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Refer to the instruction manual for steriliza-
tion instructions.

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Do not store liquids on or above the unit.

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1 Important information for users of KARL STORZ devices	V
2 Photographs of the equipment	VI
3 Controls, displays, connectors, and their uses	VIII
4 Symbols employed	X
5 General information	4
5.1 Description of the system	4
6 Safety instructions	6
6.1 Explanation of warnings and cautions	6
6.2 Intended use	9
6.3 User qualifications	9
6.4 Safety precautions at the site of installation	10
6.5 Safety precautions when operating the system	10
7 Operating instructions	11
7.1 Unpacking	11
7.2 Basic equipment	11
7.3 Installing and connecting the unit	12
7.4 KARL STORZ SCB	13
7.5 Connecting the footswitch	13
7.6 Connecting application instruments (Motor or shaver handpiece)	14
7.7 Coolant and irrigation pump	14
7.8 Preparing the tubing set	15
7.9 Insert the pump tubing	16
7.10 Initial operation	20
7.11 UNIDRIVE® S III function	20
7.12 Selecting the language	21
7.13 Selecting the motor connection	22
7.14 Overview of application instruments	23

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8 Configuring the operating mode.....	24
8.1 Select the 'EC Mikro Motor' operating mode..	24
8.1.1 Application handpiece: drill.....	24
8.1.2 Application handpiece: intranasal drill.....	28
8.1.3 Application handpiece: Universal saw	29
8.1.4 Application handpiece: Dermatomes.....	34
8.2 Operating mode: DrillCut-X	35
8.2.1 Application handpiece: Sinus drill	35
8.2.2 Application handpiece: Shaver	39
8.3 Operating mode: DrillCut-X II	40
8.3.1 Application handpiece: Sinus drill	40
8.3.2 Application handpiece: Shaver	44
8.3.3 Application handpiece: Shaver Plus	45
8.4 Using the second motor connection	46
8.5 Mode: Settings	47
8.5.1 System information	47
8.5.2 Error messages	48
9 Maintenance	49
9.1 Changing the fuse	49
9.2 Cleaning, disinfection and sterilization	50
9.3 Cleaning, disinfection and sterilization of reusable tubing.....	51
9.3.1 Preliminary cleaning of reusable tubing	52
9.3.2 Preparation	53
9.4 Sterilization	55
9.4.1 Steam sterilization.....	55
9.5 Fractionated prevacuum procedure	56
9.6 Maintenance and safety check	57
9.6.1 Maintenance.....	57
9.6.2 Safety check.....	57
9.7 Service and repair.....	59
9.8 Disposal.....	59
9.9 Repair program.....	60
9.10 Important information.....	60

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9.11 Limitation of liability.....61	
9.12 Manufacturer's warranty.....61	
10 Technical description.....62	
10.1 General.....62	
10.2 Error table UNIDRIVE®.....64	
10.3 Troubleshooting.....66	
10.4 Technical Data.....68	
10.5 Standard conformity.....70	
10.6 Directive compliance.....70	
10.7 Technical documentation.....71	
10.8 Block diagram.....72	
11 Spare parts, recommended accessories 73	
11.1 Spare parts.....73	
11.2 Recommended Accessories.....73	
12 Appendix.....76	
12.1 Cleaning agents and disinfectants.....76	
12.1.1 Manual preparation.....76	
12.2 Electromagnetic Compatibility (EMC) Information.....80	
13 Warranty cards.....94	
14 Subsidiaries.....93	

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5 General information

5.1 Description of the system

The new KARL STORZ motor system UNIDRIVE® S III ENT is characterized by the most modern electronics, flexibility and efficiency.

The system can be applied universally with brushless high-powered electronic drives integrated in the instrument or separately in the S III ENT range of application.

All settings and functions can be controlled via the integrated touch-screen. The settings and functions appear as "Text" in the display.

The defaults are stored, i.e. they only have to be selected one time (Exception: clockwise/counter-clockwise rotation; the system always starts in the clockwise setting.)

By simply activating the respective menu field, e.g. device information can be displayed. Troubleshooting in OR is easier since any faults that may occur, such as a motor fault, are displayed automatically.

Two connection sockets allow two motors to be simultaneously connected, for example, two drill handpieces in parallel.

Special features

- Automatic detection of connected motors
- Stageless variable speed control
- 2 motor outputs
- Max. speed can be limited as default setting
- Constant high motor power throughout the total speed range
- Processor controlled speed and motor torque
- Integrated coolant and irrigation pump with homogenous irrigation efficiency over the entire irrigation power range
- Optimized operator control with the touch screen
- Simple and clearly readable operating elements on a color display
- Pump activated via footswitch

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With KARL STORZ SCB connections for the KARL STORZ Communication Bus (KARL STORZ SCB)

Can be used with:

- INTRA drill handpiece (2524xx, 2525xx)
- Stammberger-CASTELNUOVO DrillCut-X shaver handpiece (40711050) with shaver attachment (412xx, 413xx)
- Stammberger-CASTELNUOVO DrillCut-X shaver handpiece (40711050) with sinus drill attachment (413xx)
- Drill Cut-X II (40712050/40712055) shaver handpiece with shaver attachment (412xx, 413xx)
- Drill Cut-X II (40712050/40712055) shaver handpiece with sinus drill attachment (413xx)
- STAMMBERGER-SACHSE Intranasal Drill (660000)
- Microsaw (254000)
- Micro sagittal saw (254100)
- Osseoscalpel (254200)
- Micro compass saw (254300)
- Dermatom (253xxx)

i NOTE: Please read the provided instructions carefully before operating the respective application instruments.

Abbreviations used:
*S III ENT = Ear Nose Throat

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Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

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6.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved.

CAUTION: A Caution indicates that particular maintenance procedures or safety precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information

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WARNING: Read this instruction manual thoroughly before operating the system. Read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

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 **WARNING:** The UNIDRIVE® S III ENT must only be used in medical facilities that conform to applicable national, state and local electrical codes.

 **WARNING:** Keep out of reach of patients.

 **WARNING:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

 **WARNING:** Combinations of medical devices are only safe if:

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

 **WARNING:** Connect the device to a properly grounded 'Hospital grade' or 'Hospital use' electrical outlet. Check grounds regularly. Routinely inspect the power cord and plug and do not use if the plug or the cord is damaged.

 **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

 **WARNING:** Do not use in the presence of flammable anesthetics. There is a risk of explosion.

 **WARNING:** The UNIDRIVE® S III ENT, in conjunction with the application instruments to be used, is only designed for intermittent operation (50%, t=1 min). (See Technical Data on p. 68).

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WARNING: Always unplug the unit before performing any maintenance, including cleaning.

WARNING: To avoid the risk of electrical shock, do not open the device. Refer servicing to the manufacturer or qualified personnel that are authorized by the manufacturer. Removal of the cover by non-authorized personnel will void the warranty.

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CAUTION: Only operate the unit with the voltage stated on the unit's identification plate.

CAUTION: When replacing fuses, use only fuses of the same rating.

CAUTION: The UNIDRIVE® S III ENT may only be used with accessories, and disposable items which have been designated for use by KARL STORZ.

CAUTION: Do not store liquids on or above the device. Do not allow fluids to get into the unit.

CAUTION: If the motor becomes blocked for any length of time, it can lead to overheating and causing damage to the motor. If the blockage cannot be disengaged using oscillating operation (if possible), the motor must be shut off and the blockage must be removed manually.

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NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

NOTE: For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices. Current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.

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Before using the UNIDRIVE® S III ENT motor system on a patient, it is imperative that the instructions in this manual are understood and are followed carefully.

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6.2 Intended use

The UNIDRIVE® S III ENT motor system is universally applicable in conjunction with the respective, suitable application instruments in the S III ENT application range.

The exact 'Intended use' of the application instruments is detailed in the respective instruction manual provided with the application instrument.

The UNIDRIVE® S III ENT motor system may only be used with accessories, parts and disposable items that have been designated by KARL STORZ as suitable for the unit.

Unauthorized conversions or modifications to the device are not allowed for safety reasons.

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6.3 User qualifications

The UNIDRIVE® S III ENT motor system in conjunction with an application instrument may only be used by physicians and/or medical assistants who have been instructed in the use of this device/instrument.

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
6.4 Safety precautions at the site of installation

The UNIDRIVE® S III ENT motor system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The UNIDRIVE® S III ENT is equipped with a connector for the ground. Connect this in accordance with the state, local and national electrical codes.

The UNIDRIVE® S III ENT motor system is not intended for use in hazardous areas. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

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6.5 Safety precautions when operating the system

! WARNING: The system must be tested for operational safety and proper condition prior to each surgical procedure.

During surgery with the UNIDRIVE® S III ENT, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

CE marked equipment has been tested and is found to comply with the EMC limits of EN/IEC 60601-1-2 :2001 CISPR 11 Class B.

Please read the Electromagnetic Compatibility Information in the appendix (p. 80-92).

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7 Operating instructions

7.1 Unpacking

Carefully remove the unit and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can be used if the unit has to be transported.

7.2 Basic equipment

KARL STORZ Set 407016 01-1

- 1 UNIDRIVE® S III ENT 407016 01-1
- 1 Power cord 400 A
- 1 Two-pedal footswitch with proportional function (200166 30)
- 1 Tubing set for irrigation, can be sterilized, 207116 40
- 1 Clip-Set 207116 21
- 1 3 pcs. mtp Tubing Sets, single use, sterile 031131-01
- 1 Irrigator stand 207116 80
- 1 SCB connecting cable, length 100 cm, 200901 70
- 1 Instruction manual 96056025D

i NOTE: Optional accessories are described in section 'Recommended accessories'.

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7.3 Installing and connecting the unit


NOTE: The UNIDRIVE® S III ENT including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

! WARNING: The UNIDRIVE® S III ENT is not intended for use in hazardous areas. Do not operate the system within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
Place the UNIDRIVE® S III ENT on a level surface. The UNIDRIVE® S III ENT is equipped with a connector ⑨ for attaching a ground wire. The unit's ground should be installed by a qualified electrician.

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Before connecting the UNIDRIVE® S III ENT to the power supply line, verify that the line voltage stated on the unit's identification plate conforms with that of the power supply lines to be used.


Connect the power, insert the power cord end into the receptacle ⑪ of the unit.

! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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
7.4 KARL STORZ SCB

① NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector is equipped with a protection mechanism. Pull the protection mechanism of the SCB connector and insert the plug into one of the SCB sockets ⑫. Connect the other end of the cable to a KARL STORZ SCB control unit (KARL STORZ Communication Bus) or other SCB devices (see instruction manual KARL STORZ SCB System).

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7.5 Connecting the footswitch

To use the footswitch ⑬, insert the connecting cord* of the footswitch into the socket for the footswitch ⑦ on the back of the UNIDRIVE® S III ENT. Align the red dot on the connecting cord with the red dot on the footswitch socket.

① NOTE: If no footswitch is connected when the unit is switched on, the message: 'E11: NO FOOTSWITCH' flashes in the status line of the display (see arrow).

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
7.6 Connecting application instruments (Motor or shaver handpiece)

Prepare the application instrument according to the instruction manual.

NOTE: Please read the instruction manual provided with the application instrument carefully!

The connecting cable* coming from the application instrument (EC Micromotor or shaver handpiece) must be connected to motor output M1 ③ and/or M2 ④.

The connection will be made properly when the red dot on the connection cable is aligned with that on the socket.

NOTE: If no 'motor' is connected when the unit is switched on, the message: 'E12: NO MOTOR' flashes in the status line of the display.

* NOTE: To unplug the connecting cable, pull out the plug – do not pull directly on the cable!

7.7 Coolant and irrigation pump

The UNIDRIVE® S III ENT control unit is equipped with an integrated roller irrigation pump which takes on the function of a coolant pump as well.

NOTE: Use physiological saline solution as a coolant and irrigation fluid.

The pumping function is available with the following application instruments: SHAVER, DrillCut-X II, DRILLCUT-X, DRILL, SINUS BURR, INTRANASAL and SAWS.

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
In some cases, the irrigator stand should be mounted

If necessary, insert the irrigator stand for hanging an irrigation bag or bottle in the holder ⑥. Lock the irrigator stand in place by turning the knob on the holder (clockwise).

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CAUTION: Mount the irrigation stand in such a way that the irrigation bag neither hangs directly above the device nor touches the device (heat emission) (see section 'Safety instructions').


NOTE: If an irrigation bottle is used, please ensure that the bottle is ventilated. Otherwise irrigation performance decreases.

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7.8 Preparing the tubing set

NOTE: The tubing set is not sterile when delivered. For sterilization, see pages 55 ff.

NOTE: The tubing is delivered ready-assembled. If the tubing has been disassembled, for preparation for example, it must be reassembled as described below.

Connect both ends of the pump tubing to the connectors until they click into place.

Attach inflow tubing to the narrow pump tubing connector (large diameter). Attach outflow tubing to the wide pump tubing connector (small diameter).

After attaching the pump tubing, press the movable clamp ring on the wide pump tubing connector against the pump tubing as far as it will go.

WARNING: It is essential to carry out a leakage test every time before you attach the pump tubing to the device. See page 54.

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7.9 Insert the pump tubing

Flip the lever for retracting/extending the pump roll down into the horizontal position.

! **CAUTION:** Use only the tubing set that is intended for use with the UNIDRIVE® S III ENT.

i **NOTE:** The pump tubing connectors have different widths so that the pump tubing can only be inserted in one direction. There is also a directional arrow on the connector to the drainage tubing.

Attach the pump tubing to the roller pump and press both connectors into the fastening grooves.

i **NOTE:** If the connectors cannot be pressed into the fastening grooves, they have not been attached correctly to the tubing. If necessary, attach the connectors on the pump tubing more securely.

i **NOTE:** Due to application wear, we recommend changing the pump tubing after 5 uses.

! **CAUTION:** If the pump tubing is perforated, there is a danger of fluid entering the device!

Move the lever for retracting/extending the pump roller into the vertical position and the pump tubing is locked in place.

i **NOTE:** Without locking the pump tubing, the pump cannot be started, error message 'E16: IRRIGATION PUMP NOT CLOSED' appears.

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
Removing the tubing set

Press lever down to the horizontal position.
Remove the pump tubing from the reel pump.

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The inflow tubing of the silicone tubing set for
irrigation is equipped with a puncturing pin and a
tubing clamp.


Close the tubing clamp.

Insert a piercing punch into the intended position
on the irrigation bag or irrigation bottle.

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If necessary, attach a cooling/irrigation channel
onto the drill handpiece/microsaw (see arrow).



Attach the other end of the silicone tubing to the
irrigation connection of the DrillCut-X shaver hand-
piece, the drill handpiece or the microsaw.

Using clips from the clip set (Art. no. 20711021),
attach the silicone tube to the DrillCut-X and
DrillCut-X II drill handpiece/micro saw and/or con-
necting cable.

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**Disposable irrigation tubing set (mtp Tubing
Set 031131-10)**

Insert the tubing into the holder.
Connect the relevant tubing end with the irrigation
bottle (puncture needle) or attach to the irrigation
connector of the DrillCut-X shaver handpieces, drill
handpiece or the micro-saw.

NOTE: If an irrigation bottle is used open
the ventilation cap at the puncture needle in
order to ventilate the bottle.

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Irrigation pump control

NOTE: Operating the pump functions via the
touch screen is described in the respective
section under 'Initial operation'.

The pump function is either disabled or enabled
with the 'Pump' operating button on the touch
screen of the control unit. If the pump function is
enabled (pump button turns orange), the pump is
automatically activated with the motor when the
right footswitch pedal is activated and it is ...

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... deactivated when the pedal is released.

NOTE: Activating the left-hand footswitch
pedal for a longer period of time activates
the irrigation pump alone (if it was activated
previously).

The pump runs backwards for a short time after
deactivating to aspirate the irrigation solution. This
prevents any liquid from dripping.

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Setting the flowrate – Touch screen

The '+/-' keys in the touch screen's 'Pump' screen can be used to set the specified irrigation flow rate for the pump in the respective operating mode in 9 pump levels.

The currently selected pump level is displayed accordingly.

NOTE: If the pump level is set to '0', the pump does not start.

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Setting the flowrate – Footswitch

NOTE: The irrigation flowrate can also be set using the footswitch.

Pressing the left footswitch pedal briefly raises or lowers (depending on setting) the irrigation flow rate level by level.

NOTE: Setting with the left foot pedal is only possible if the right foot pedal is not active; i.e. not while drilling or while operating the motor/shaver.

NOTE: The last value that was set for the solution flowrate while in the defined operating mode is retained even after switching the control unit off.

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7.10 Initial operation

Switch the UNIDRIVE® S III ENT on at the power supply ① (power switch is illuminated, screen turns white for a moment).

- ❗ **NOTE:** Before switching the unit on, the user must connect the footswitch and the appropriate handpiece. The software is designed to verify the connections and will indicate an error message in the lower screen margin, if necessary.

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The screen shown on the left appears for the time it takes the application to start.

- ❗ **NOTE:** The next pages contain a description of the simplest application example; adopting the most recently set mode, as well as the adjustments for speed and pump settings. During the switch-on and switch-off procedures, a single signal tone is heard..

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7.11 UNIDRIVE® S III function

A start screen is displayed once the unit has started. The function only becomes active once application instruments have been connected with the unit.

- ❗ **CAUTION:** To prevent damage, do not use any object other than fingertips to activate the sensitive touch screen.

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The screen displayed shows the start screen with EC micromotor connected.

- ❗ **NOTE:** Error messages are shown in text in the status line of the operation screen if necessary (flashing text).

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7.12 Selecting the language


It may be required to define a different language for the operation screen since the device is set to English as the default in the factory.

This is done by pressing the 'menu' symbol (wrench) in the bottom left.

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Activate the 'LANGUAGE' button (2nd button from the top).

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X = Aktion abbrechen

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Select the desired language*.

Select one of the ten available languages by pressing the respective language button. The languages currently available are displayed over two screens.

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7.13 Selecting the motor connection

Various application instruments operated within various speed ranges can be connected up to the control unit. To do so, the required instrument can be optionally connected to the motor connection socket (M1 and M2). Having detected the handpiece automatically, the UNIDRIVE® S III loads the standard parameters and displays them on the screen. The symbol for the respective instrument is displayed in the bottom right of the screen.

There are two connection sockets available, which can be used with different application instruments.

The motor connection socket can be selected via the buttons for the respective application in the lower right of the screen.

It is only possible to switch between the individual motor connection sockets when both channels are occupied.

ⓘ NOTE: It is also possible to switch between M1 <-> M2 by pressing the button **Ⓜ** on the footswitch for longer (> 2 s) (see also foot-switch function 1.5).

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
7.14 Overview of application instruments

1. INTRA drill handpiece (2524xx) with EC Micromotor (20711032) and drill (26xxxx, 649xxx K, 649xxx GK, 649xxx GL, 649xxx L)	Page 24	DRILL	1,000 – 60,000 rpm
2. Drill Cut-X II shaver handpiece (40712050/40711050) with shaver attachment (412xx, 413xx)	Page 40	DrillCUT-X II	300 – 10,000 rpm Preset = 4,000 rpm
3. Drill Cut-X II shaver handpiece with sinus (40712050/40711050) (413xx) drill attachment	Page 40	DrillCUT-X II	300 – 12,000 rpm Preset = 12,000 rpm (max.)
4. Stammberger-CASTELNUOVO DrillCut-X sh. handp. (40711040) with shaver attachment (412xx, 413xx)	Page 35	DRILLCUT-X	300 – 7,000 rpm Default = 7,000 rpm (max.)
5. Stammberger-CASTELNUOVO DrillCut-X shaver handpiece (40711040) with Sinus burr attachment (413xx)	Page 35	SINUS BURR	800 – 12,000 rpm (40711040) Default = 12,000 rpm (max.)
6. Microsaw (254000)*, Micro-Sagittalsaw (254100)*, Osseoscalpel (254200)*, Micro-piercing saw (254300)* * with EC Micromotor (20711032)	Page 45	SAWS	1,000 – 20,000 rpm Default = 15,000 rpm
7. Dermatomy (253xxx) with EC Micromotor (20711032)	Page 34	DERMATOME	1,000 – 8,000 rpm

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
8 Configuring the operating mode

After switching on and enabling the functions of the UNIDRIVE® S III ENT, the adjacent screen appears.

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8.1 Select the 'EC Mikro Motor' operating mode

(INTRA drill handpiece 2524xx with EC micromotor 20711032)

To activate the 'EC micromotor' operating mode, connect the EC micromotor up to the UNIDRIVE® S III ENT using the connection cable.

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
Selecting the 'Drill' application.

All applications operated with the EC micromotor are available for selection here. These are: drill, intranasal drill, micro compass saw, and dermatome.

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
Press the confirm button in the bottom right to be able to use the selected function.

NOTE: Once an application has been selected, the display can be returned to the run screen using the confirm button.

The term 'Drill' can now be seen in the status display (top left).

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Configuring the operating mode: Activating the 'SOFT START' option

The option 'SOFT START' makes it possible to accelerate the drill handpiece continuously to the default maximum speed when activating the foot-switch pedal within a predefined amount of time.

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The left-hand field on the operation screen is used for setting and displaying the motor rpm. The available speed range is shown in the left-hand bar indicator at the top and bottom (in this case 1,000 to 40,000 rpm).

Use the '+/-' buttons to define the maximum speed (set value). This preset speed is displayed by the orange bar and figures.

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NOTE: Please observe the maximum permitted rpm for the selected attachment!

If the EC micromotor is switched on via the right footswitch, the bargraph display for the speed turns light blue and displays the current speed (actual value).

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Configuring the operating mode: Configuring the coolant and irrigation pump

In every operating mode (with the exception of 'DERMATOMES'), the integrated cooling and irrigation solution pump can be activated as an option. To activate, press the 'PUMP' button in the center of the screen.

NOTE: Information on adjusting the irrigation pump settings can be found on page 19.

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If the pump is activated, the button and the bargraph display are also activated.


The irrigation flow rate can be preset using the '+/-' keys (setpoint value).

The setpoint value is displayed by the blue line in the bargraph display.

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If the handpiece is started with the right-hand footswitch, the pump starts automatically and supplies the defined flow quantity.

The left-hand footswitch pedal (actuation >2 s) only starts the pump and the irrigation solution is supplied with the defined irrigation flow. The handpiece is not activated.

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The sub-menu allows you to choose from the saws on offer.

The following saws are available to you:

- Micro compass saw (as of page XX)
- Oscillating micro compass saw (as of page XX)
- Osseoscalpel (as of page XX)
- Sagittal saw (as of page XX)

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
**Application handpiece:
Micro compass saw**


You are now in the 'Universal saw' sub-menu.

NOTE: Information on how to access this sub-menu can be found on p. 29.

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Select the 'Micro compass saw' application.

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The term 'Micro compass saw' can now be seen in the status display (top left).

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
Once the confirm button has been pressed, the 'Micro compass saw' function is available. The correct settings for this function are loaded automatically.

To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 24.

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**Application handpiece:
Oscillating micro compass saw**
You are now in the 'Universal saw' sub-menu.

NOTE: Information on how to access this sub-menu can be found on p. 29.

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
Select the 'Oscillating micro compass saw' application.

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The term 'Oscillating micro compass saw' can now be seen in the status display (top left).

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Once the confirm button has been pressed, the 'Oscillating micro compass saw' function is available. The correct settings for this function are loaded automatically.

To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 24.

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**Application handpiece:
Osseoscalpel**


You are now in the 'Universal saw' sub-menu.

NOTE: Information on how to access this sub-menu can be found on p. 29.

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
The term 'Osseoscalpel' can now be seen in the status display (top left).

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Once the confirm button has been pressed, the 'Osseoscalpel' function is available. The correct settings for this function are loaded automatically. To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 24.

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**Application handpiece:
Sagittal Saw**


You are now in the 'Universal saw' sub-menu.

NOTE: Information on how to access this sub-menu can be found on p. 29.

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Select the 'Sagittal saw' application.

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The term 'Sagittal saw' can now be seen in the status display (top left).

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Once the confirm button has been pressed, the 'Sagittal saw' function is available. The correct settings for this function are loaded automatically.

To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 24.

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8.1.4 Application handpiece: Dermatomes

To change the operating mode (selection of application handpiece), press the button of the displayed attachment on the right of the screen.

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
Select the 'Dermatomes' application.

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The term 'Dermatomes' can now be seen in the status display (top left).

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
Once the confirm button has been pressed, the 'Dermatomes' function is available. The correct settings for this function are loaded automatically.

To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 24.

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8.2 Operating mode: DrillCut-X

To activate the 'DrillCut-X' operating mode, connect the DrillCut-X up to the UNIDRIVE® S III ENT using the connection cable.

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
8.2.1 Application handpiece: Sinus drill

In order to be able to work with the selected application handpiece, please insert the unit into the selected motor connection socket. The handpiece is detected automatically and the corresponding screen loaded.

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Select the 'Sinus drill' application.

All applications operated with the DrillCut-X are available for selection here.

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The term 'Sinus drill' can now be seen in the status display (top left).

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The top field of the run screen is used for setting and displaying the motor speed. The available speed range is shown in the left-hand bargraph display at the top and bottom (in this case 300 to 10,000 rpm). The maximum speed (setpoint value) can be preset using the '+/-' buttons. This preset speed is displayed by the orange bar and figures in the bargraph display.

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NOTE: The 'start-up value' for the max. speed is 3,000 rpm or less.

If the 'handpiece' is switched on via the right foot-switch pedal, the bargraph display and the figures for the speed turn light blue and display the current speed (actual value).

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The integrated coolant and irrigation pump can be optionally activated in all operating modes (apart from 'DERMATOMES'). For the activation, the 'PUMP' button on the right of the run screen must be clicked.

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The irrigation flow can be preset using the '+/-' keys' (setpoint value).

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Footswitch function

The footswitch has the following functions in the 'SHAVER' operating mode:

a) Activate the left pedal

only the pump is activated (irrigation), flow according to defined value (see Activating the pump '2.4')

b) Short, quick actuation of the left pedal

The irrigation flow rate can be increased or reduced with the left pedal as of V1.03 (depending on the setting in the 'Accessories configuration' menu).

c) Activate the right pedal

this pedal has a proportional function (Gas pedal function), the speed is varied depending on the pedal position. The handpiece is started, pump runs with the defined irrigation flow quantity.

d) Press button ⑭ for longer than 2 sec.

Switches from motor output 1 (M1) and motor output 2 (M2), or vice versa.

e) Press left and right pedal simultaneously

Pushing both pedals simultaneously allows you to switch between the shaver and sinus drill modes. A signal tone is emitted when you switch between modes.

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The term 'Shaver' can now be seen in the status display (top left).

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
Once the confirm button has been pressed, the 'Shaver' function is available. The correct settings for this function are loaded automatically.

To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 38.

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8.3 Operating mode: DrillCut-X II

To activate the 'DrillCut-X II' operating mode, connect the DrillCut-X II up to the UNIDRIVE® S III ENT using the connection cable.

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8.3.1 Application handpiece: Sinus drill

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
Selecting the 'Sinus drill' application.
All applications operated with the DrillCut-X are available for selection here.
These are: sinus drill, shaver.

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NOTE: Once an application has been selected, the display can be returned to the run screen using the confirm button.
The term 'Sinus drill' can now be seen in the status display (top left).

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The option 'SOFT START' makes it possible to accelerate the drill handpiece continuously to the default maximum speed when activating the foot-switch pedal within a predefined amount of time. Press the 'SOFT START' button in the top right of the screen.

The chosen settings will be saved when the required time is pressed. You are then returned directly to the application screen.

NOTE: The selected 'ramp' (start-up behavior) is used for the respective function. It is also generally possible to use different ramps for different applications.

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Configuring the operating mode: Preselecting the maximum speed (setpoint value)

The top field of the run screen is used for setting and displaying the motor speed. The available speed range is shown in the left-hand bar indicator at the top and bottom (in this case 300 to 12,000 rpm). The maximum speed (setpoint value) can be preset using the '+/-' buttons. This preset speed is displayed by the orange bar and figures in the bargraph display.

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
8.3.3 Application handpiece: Shaver Plus

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Select the 'Shaver Plus' application.

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И **NOTE:** Once an application has been selected, the display can be returned to the run screen using the confirm button.

The term 'Shaver Plus' can now be seen in the status display (top left).

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Once the confirm button has been pressed, the 'Shaver Plus' function is available. The correct settings for this function are loaded automatically.

To adjust the individual settings in this mode, please proceed as described in the 'Configuring the operating mode' section as of page 43.

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8.4 Using the second motor connection

Since the UNIDRIVE® S III ENT has two motor connections (channels), a second application handpiece can be connected.

The configurations and settings described in this installation manual apply for both connections (M1 and M2). The two motors can be configured independently of one another.

To switch between the motor output 1 (M1) and motor output 2 (M2), the corresponding button can be selected in the display (bottom right) or the button ⑭ on the footswitch must be pressed for longer than 2 sec.

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8.5 Mode: Settings


8.5.1 System information

To call up the menu functions, press the 'Menu' button (wrench) in the bottom right.

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
This opens the 'MENU' screen (image on the left). Press the 'SYSTEM INFORMATION' button to obtain detailed information on your UNIDRIVE® S III ENT.

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The device number (serial number) is displayed at the top.


The hours of operation and the software version is shown in the lower part of the screen.

Activating the 'OPERATING MODE' button returns you to the operation screen, from which the MENU was called. This applies for every menu level.

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8.5.2 Error messages

By pressing the 'SERVICE' button from the 'MENU', the service sub-menu is opened.

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In the top portion of the screen, the device number (serial number) is displayed. In the lower portion of the screen, the error messages that have been shown are displayed in the order of importance.

E01-E10: Hardware error messages, which could possibly lead to a defective device.

E11-E30: User error messages that lead back to an incorrect configuration of the device (e.g. footswitch not connected). See Error table on p. 64.

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Three error messages are shown in the example on the left:

E 11: NO FOOTSWITCH (no footswitch connected)

E 12: NO MOTOR IN M1 (no motor/handpiece connected to M1)

E 16: IRRIGATION PUMP NOT CLOSED (the pump head is not closed)

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9.2 Cleaning, disinfection and sterilization

-  **NOTE:** Please read the instruction manual provided with the application instrument carefully!
-  **NOTE:** The instructions for 'cleaning, sterilization and care of KARL STORZ instruments' must be followed. The procedures for cleaning, disinfecting and sterilization are explained in detail there.
-  **NOTE:** For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices. The current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.

 **WARNING:** Before any cleaning work is performed, the device is to be disconnected from the power source.

 **CAUTION:** Make absolutely sure that no liquid can get into the device housing.

Wipe-clean the exterior surfaces of the unit with a disposable cloth moistened with disinfectant. For wipe-down disinfection, we recommend using only surface disinfectants (no solutions predominantly containing alcohol such as rapid-action disinfectants) and 'disposable disinfectant methods' based on ammonium chloride.

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9.3 Cleaning, disinfection and sterilization of reusable tubing

! WARNING: Reusable tubing must be cleaned, disinfected and sterilized before use and subsequent reuse using validated preparation procedures.

! WARNING: When cleaning and disinfecting contaminated instruments (tubing sets), the local guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed.

! WARNING: Prior to the sterilization of reusable tubing, particularly in cases of manual preparation, it must be ensured that all chemical residue has been rinsed off or removed (see also the preparation steps).

! WARNING: Disposable tubing sets should not be reused.

i NOTE: For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices. The current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.

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9.3.1 Preliminary cleaning of reusable tubing

Place the reusable tubing in a container filled with a cleaning solution (as per the manufacturer's instructions) immediately after use to prevent contaminants from drying on the surface of the tubing set.

Water quality requirements

Tap water can be used to prepare the cleaning and disinfectant solutions.

When rinsing for the last time, for example, to remove any chemical residue, microbiologically pure/sterile water must be used.

CAUTION: When preparing and using the solutions, follow the manufacturer's instructions paying close attention to proper concentration and exposure times. Prolonged immersion can lead to material changes.

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
9.3.2 Preparation

Manual preparation of reusable tubing

WARNING: For all cleaning and disinfection procedures on contaminated instruments, the guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed.

1. Disconnect the reusable tubing and immerse it in cleaning solution (as per the manufacturer's instructions).
2. Rinse the tubing through under the (water) surface using a syringe. Repeat the procedure several times if necessary. Remove soiling with brushes, sponges, and lint-free disposable cloths. Cleaning accessories are available from KARL STORZ (see catalog).
3. After cleaning, rinse the tubing and place it in a suitable container with disinfectant solution, making sure that all air bubbles have been expelled. Ensure that the tubing is completely covered in the solution for the duration of exposure.
4. Cleaning brushes must be cleaned and disinfected immediately.
5. Remove the tubing after the exposure time and rinse with microbiologically pure/sterile water. Ensure that all chemical residue in and on the tubing is removed.
6. Dry the tubing with sterile compressed air.
7. After preparation, inspect the tubes for cleanliness and damage. If contaminants can still be detected, the preparation procedure must be repeated.

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Machine preparation of reusable tubing


To prepare using a machine, connect the reusable tubing to the retaining mechanisms in such a way that the instrument channels can be completely rinsed through.

Visible, heavy soiling must be removed in preliminary cleaning prior to machine preparation.

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
Assemble tubing set as described in the section 'Preparing the tubing set' on p. 15, making sure that connections are secure to prevent leakage in the connector area.

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! WARNING: Check the entire tubing set for leaks after preparation. As a result, it is possible to detect a leak early, through which liquid could escape and enter the unit.

To check for leaks use leakage tester 13242 XL for example. For this purpose connect one end to the leakage tester and block off the other end, with your finger for example. If the tubing set develops a leak, a pressure drop will be observed on the pressure gauge of the leakage tester.

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9.4 Sterilization

We recommend steam sterilizing the reusable tubing before use at 134°C (+3°C) using the fractionated prevacuum procedure.

WARNING: Successful sterilization is only possible on clean, disinfected surfaces.

WARNING: The required sterilization parameters only apply in conjunction with a properly maintained and validated sterilization apparatus.

WARNING: Deviations from the recommended sterilization parameters are to be confirmed by the user.

9.4.1 Steam sterilization

KARL STORZ has validated the following sterilization method: Fractionated prevacuum method (steam sterilization)

NOTE: Use the fractionated vacuum method (acc. to DIN EN 550/ISO 11134), taking into account the respective applicable national requirements.

The containers are to be placed in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

Load the sterilizer according to its instruction manual.

At the completion of the steam sterilization cycle, all instruments should remain untouched until adequately cooled.

WARNING: Burn hazard! Parts are hot after steam sterilization. Allow them to cool down.

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9.5 Fractionated prevacuum procedure

The fractionated prevacuum procedure also consists of the four phases mentioned previously.

- **Conditioning phase**

A vacuum is created up to four times in the conditioning phase. Then, the steam is blown in.

- **Sterilization phase**

The sterilization is performed at a temperature of 134°C (+3°C), a pressure (absolute) of 300 kPa (3 bar) for a duration of 5 minutes (minimum time).

- **Evacuation phase**

The steam is released.

- **Drying phase**

Drying takes place under renewed application of a vacuum over a period of approx. 5–20 minutes.

Validated parameters for the fractionated prevacuum method:	
Temperature:	134°C (+3°C)
Pressure (p _{absolute}):	300 kPa (3 bar)
Exposure time:	5 minutes for all thermoresistant medical devices

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9.6 Maintenance and safety check**9.6.1 Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

9.6.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once per year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

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Electric measurements:

- Inspecting the device safety fuses
- Protective ground resistance measured according to IEC 62353 ed1.0 Limit value: 0.2 ohms
- Earth leakage current measured according to IEC 62353 ed1.0: Limit value: NC 500 µA, SFC: 1000 µA
- Patient leakage current measured according to IEC 62353 ed1.0: Limit value: 100 µA (BF), 10 µA (CF)
- Leakage currents must not exceed the limit values specified above.

Test for proper functioning:

A functional check must be carried out in accordance with the instruction manual.

The safety check must be entered in the unit log-book and the results of the check must be documented.

If the device does not function reliably and/or safely, it must be repaired.

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9.7 Service and repair

Repairing defective devices is only to be performed by personnel that have been authorized by KARL STORZ and using original KARL STORZ parts.

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9.8 Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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9.9 Repair program

Usually to bridge the repair period, you will receive a unit on loan, which you will then return to KARL STORZ as soon as you receive the repaired unit.

In Germany, you can refer the repairs directly to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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9.10 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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9.11 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

9.12 Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the warranty card and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

i NOTE: The unit must always be sent to the following address (even during the guarantee period; where applicable, with guarantee card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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10 Technical description

10.1 General

In the event of the unit malfunctioning, the troubleshooting list provided in this chapter will assist the user in finding the possible cause. The troubleshooting list indicates how to find and rectify the various malfunctions.

WARNING: High voltage is present in the unit's interior. The unit must only be opened by an authorized KARL STORZ customer service technician.

WARNING: Maintenance work that is not mentioned in this chapter must only be performed by technical personnel that are authorized by KARL STORZ.

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General information

Error messages are output directly to the screen on the UNIDRIVE® S III ENT. Send the device in for servicing when a hardware error message is displayed (E01-E10). Take appropriate action for user error messages (E11-E30). E.g. E11 'NO FOOTSWITCH' – check whether the footswitch is properly connected. Detailed information can be obtained from the error table on page 64 (see section 'Service mode' as well).

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10.2 Error table UNIDRIVE®

Error message	Cause	Remedy
1 = E01: MEMORY ERROR	Parameters could not be saved successfully.	Turn device on/off. Load default values. Or send unit in for repair.
2 = E02: MOTORBOARD DEFECTIVE	Motor control overloaded.	Turn device on/off. Or send unit in for repair.
3 = E03: HARDWARE CONFIGURATION ERROR	Hardware configuration is incorrect.	Send unit in for repair.
4 = E04: SOFTWARE CONFIGURATION ERROR	Software configuration is incorrect.	Send unit in for repair.
5 = E05: -		
6 = E06: -		
7 = E07: -		
8 = E08: -		
9 = E09: MOTOR DEFECTIVE M1	Motor on channel 1 is defective.	Release footswitch and restart. If the error persists, a different handpiece must be used.
10 = E10: MOTOR DEFECTIVE M2	Motor on channel 2 is defective.	See E09
11 = E11: NO FOOTSWITCH	Footswitch not connected correctly.	Connect footswitch. If the error persists, check whether the error is on the footswitch or unit. Send the defective component in for repair.
12 = E12: NO MOTOR	No motor/handpiece connected	Connect motor/handpiece.
13 = E13: TOO HIGH TEMPERATURE MOTOR M1	Motor/handpiece on channel 1 overheated.	Allow the connected handpiece to cool down or operate an additional handpiece on the 2nd channel.
14 = E14: -		
15 = E15: TOO HIGH TEMPERATURE MOTOR M2	Motor/handpiece on channel 2 overheated.	See E13
16 = E16: IRRIGATION PUMP NOT CLOSED	'Pump lever' not connected correctly.	Set pump lever vertical. If the error persists, send unit in for repair.
17 = E17: MOTOR NOT SUPPORTED M1	The motor/handpiece connected on channel 1 is not supported.	Only use handpieces approved by KARL STORZ. A service update may be necessary.
18 = E18: MOTOR NOT SUPPORTED M2	The motor/handpiece connected on channel 2 is not supported.	See E17
19 = E19: MOTOR STOPPED M1	Motor/handpiece on channel 1 blocked in operation.	Release footswitch and restart. If the error persists, clean the handpiece or exchange the attachment.
20 = E20: MOTOR STOPPED M2	Motor/handpiece on channel 2 blocked in operation.	See E19
21 = E21: NO SCB MODULE	Unit is not connected to an active SCB bus	Insert SCB cable and activate SCB bus.

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Power switch lamp will not light up.

Possible cause

- 1 Power failure.
- 2 Defective power supply fuse.
- 3 Power switch lamp defective.

Recommended remedy

- 1 Check that there is electricity to the wall outlet.
- 2 Change fuses as described in the instruction manual. Make sure you fit the correct type of fuse.
- 3 Notify customer service.

Problem

- 1 Power switch lamp is illuminated.
- 2 No indication of rotation.
- 3 No energy output when operating the foot-switch/buttons.
- 4 Speed cannot be regulated.

Possible cause

- 1 Motor control electronics faulty.
- 2 Connection to motor interrupted.
- 3a Faulty handpiece.
- 3b Handpiece blocked or jammed.

Recommended remedy

- 1, 2, 4 Notify customer service.
- 3a Connect another handpiece.
- 3b Clean the handpiece.

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Problem

Tubing set connectors do not fit the pump.

Possible cause

- 1 Tubing set is assembled incorrectly and upper and lower connectors are the opposite way around.
- 2 Connectors are not attached to the tubings properly.

Recommended remedy

- 1 Reassemble tubing set. (Upper and lower connectors are different sizes.)
- 2 Press connectors firmly in the direction of the pump tubing.

Problem

Irrigation fluid is not flowing.

Possible cause

- 1 Irrigation liquid bag empty.
- 2 Tubing clamp impedes flow.
- 3 Kinked tubing.
- 4 Irrigation bottle not aerated.
- 5 Pump is deactivated.

Recommended remedy

- 1 Use new irrigation liquid bag.
- 2 Loosen tubing clamp.
- 3 Change run of tubing.
- 4 Aerate irrigation liquid bottle.
- 5 Activate pump.

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10.4 Technical Data

UNIDRIVE® S III ENT	40701620-1
Line voltage	100 – 240 V~
Power frequency	50/60 Hz
Power fuse	2 x T 4.0 AL 250 V, IEC 127 Format
Power consumption	150 – 200 W
Operating temperature	10-40°C
Dimensions	305 x 165 x 280 mm (W x H x D)
Weight	5.2 kg
Operating mode (in conjunction with application instrument)	periodic intermittent operation
Relative On-time	50 %; t = 1 min
Storage/transport conditions:	
Relative humidity (non-condensing)	20%...95%
Temperature	(-)10°C...60°C
Atmosph. pressure	500 hPa...1080 hPa
EC micromotor drill	
Minimum speed:	1,000 U/min
Maximum speed:	40,000 U/min
Start-up maximum speed:	40,000 U/min
Mode:	Rotation
DrillCut-X II – Shaver	
Minimum speed:	300 U/min
Maximum speed:	10,000 U/min
Start-up maximum speed:	4,000 U/min
Mode:	Oscillate
DrillCut-X II – Sinus drill	
Minimum speed:	300 U/min
Maximum speed:	12,000 U/min
Start-up maximum speed:	12,000 U/min
Mode:	Rotation

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DRILLCUT-X

Minimum speed:	300 rpm
Maximum speed:	7.000 rpm
Switch-on maximum speed:	7.000 rpm
Rotation mode:	Oscillation

SINUS BURR with DRILLCUT-X

Minimum speed:	800 rpm
Maximum speed:	12.000 rpm
Switch-on maximum speed:	12.000 rpm
Rotation mode:	Rotation

INTRANASAL DRILL

Minimum speed:	3.000 rpm
Maximum speed:	60.000 rpm
Switch-on maximum speed:	60.000 rpm
Rotation mode:	Rotation

SAWS

Minimum speed:	1.000 rpm
Maximum speed 1:	20.000 rpm
Maximum speed 2:	15.000 rpm
Switch-on maximum speed 1:	20.000 rpm
Switch-on maximum speed 2:	15.000 rpm
Rotation mode:	Rotation

DERMATOMES

Minimum speed:	1.000 rpm
Maximum speed:	8.000 rpm
Switch-on maximum speed:	8.000 rpm
Rotation mode:	Rotation

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(for 40701420-1)

According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electrical shock:
Protection class I
- Degree of protection against electric shocks:
Applied part type BF
- Type of moisture protection:
drip-water protection as per IPX 1

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10.6 Directive compliance

Medical product of class IIa

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

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10.7 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

ⓘ NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
**11 Spare parts, recommended
accessories**
11.1 Spare parts

Part	Order number
Double-Pedal Footswitch, two-stage, with proportional function	20011630
Silicone Tubing Set, for irrigation, reusable, sterilizable	20711640
Clip Set, for tubing set	20711621
Pump tubing, package of 25, sterilizable	20330393
Power line fuse T 4.0 AL 250 V (100...240V ~)	2027690
SCB connecting cable, length 100 cm	20090170
Irrigator rod	20711680
Mains Cord (Schuko)	400 A
Mains Cord 'Hospital Grade' (USA)	400 B
Instruction Manual	96056025D
Operator's quick reference guide	96056xxxX

11.2 Recommended Accessories

(not included in delivery)

Part	Order number
Tubing Set, for single use, for use with UNIDRIVE® S III ENT, sterile, package of 10	mtp*
Universal-Spray, package of 6 sprayers	280052

* mtp medical technical promotion gmbh,
Take-Off GewerbePark 46
78579 Neuhausen ob Eck, Germany
Phone: +49 (0)7467 945 04-0
Fax: +49 (0)7467 945 04-99

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**Spare parts, recommended
accessories**

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**Recommended accessories in combina-
tion with high-performance EC micromotor
20711032 and connection cable 20711172**

(Not included in the set)

Part	Order no.
INTRA drill handpiece, angled, 18 cm, correspond- ing drills: standard, diamond, coarse diamond	252479
Rack for 12.5 cm drills, fold- ing model	280034
INTRA drill handpiece, angled, 15 cm, correspond- ing drills: standard, diamond, coarse diamond	252477
Rack for 9.5 cm drills, fold- ing model	280033
INTRA drill handpiece, angled, 12.5 cm *	252475
INTRA drill handpiece, straight, long, 10.4 cm *	252495
INTRA drill handpiece, straight, 8.7 cm *	252490
* Corresponding drills: standard, tungsten carbide, tungsten carbide cross cut, diamond, coarse diamond	
Tray holder for drill/burr sys- tem without rack	39552 A
Tray holder for drill/burr sys- tem with rack	39552 B
Micro sagittal saw	254100
Osseoscalpel	254200
Micro compass saw	254300
Oscillating micro saw	254000
Tray holder for micro saw system	39553 A

* mtp medical technical promotion gmbh,
Take-Off GewerbePark 46
78579 Neuhausen ob Eck, Germany
Phone: +49 (0)7467 945 04-0
Fax: +49 (0)7467 945 04-99

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**Recommended accessories in combina-
tion with high-performance EC micromotor
20711032 and connection cable 20711172**

(Not included in the set)

Part	Order no.
Intranasal drill	660000
Tray holder for intranasal drill system	39551 A
Dermatome 12 mm	253000
Dermatome 25 mm	253100
Dermatome 50 mm	253200
Dermatome 75 mm	253300
Tray holder for dermatome system	39554 A

**Recommended accessories in combina-
tion with DrillCut-X 40711050, DrillCut-X II
40712050, and DrillCut-X II N 40712055**

(Not included in the set)

Shaver blades: straight, 35° curved, 65° curved non- serrated, single-serrated, double-serrated, reusable or disposable (sterile), pack of 5.	Article number, see catalog
Sinus drill: oblong, round drill, diamond head 15° curved, 40° curved, 55° curved, disposable (sterile)	Article number, see catalog

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12 Appendix**12.1 Cleaning agents and disinfectants****Chemicals for the preparation of telescopes and endoscopic instrument sets****12.1.1 Manual preparation**

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU telescopes with ocular focusing, magnifiers with focusing, motors, and their handpieces, IMPERATOR drill handpieces, prismatic light deflectors. The following instruments produced by KARL STORZ are not suitable for cleaning in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass. When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to proper dilution and recommended immersion times. Regarding the microbiological effectiveness please ask the manufacturer.

The materials approved by KARL STORZ as regards material compatibility can be found on the following page.

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Typ

- 1 = cleaner
2 = disinfectant
3 = suitable for ultrasound baths
A = not suitable for flexible endoscopes
(A) = only for new generation flexible endoscopes (year of manufacture >2005)

NOTE: You will also find the current list of approved materials on the Internet at www.karlstorz.com.

NOTE: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

CAUTION: The use of CAUSTIC SODA can cause changes to the surface of aluminum parts (including coated parts), plastics, and soldered connections, and may reduce the service life of the instrument set.

For flexible endoscopes, only immerse the examination sheath, not the housing/handle part.

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12.1.2 Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a washer and disinfectant: HAMOU telescopes with ocular focusing, magnifiers with focusing, motors, and their handpieces, IMPERATOR drill handpieces, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special procedures must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable.

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Type

- 1 = cleaner
 2 = disinfectant
 3 = suitable for ultrasound baths
 A = not suitable for flexible endoscopes
 (A) = only for new generation flexible endoscopes (year of manufacture >2005)

ⓘ NOTE: You will also find the current list of approved materials on the Internet at www.karlstorz.com.

ⓘ NOTE: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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12.2 Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The UNIDRIVE® S III ENT model 40701620-1 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The UNIDRIVE® S III ENT model 40701620-1 is a Group 1 unit (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The UNIDRIVE® S III ENT model 40701620-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the UNIDRIVE® S III ENT model 40701620-1. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 201

Guidance and manufacturer's declaration – electromagnetic emissions

UNIDRIVE® S III ENT model 40701620-1 is intended for use in the electromagnetic environment specified below.
 The customer or user of the UNIDRIVE® S III ENT model 40701620-1 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group1	UNIDRIVE® S III ENT model 40701620-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	The UNIDRIVE® S III ENT model 40701620-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions acc. to IEC 61000-3-2	Class D (see NOTE below)	
Voltage fluctuations/ flicker emissions acc. to IEC 61000-3-3	not applied (see NOTE below)	
 NOTE: The measurement of harmonic emissions is not required (power < 75 W) The measurement of voltage fluctuations/flicker is not required (activated manually, unlikely to produce significant voltage fluctuations or flicker)		

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Table 202

Guidance and manufacturer's declaration – electromagnetic immunity

UNIDRIVE® S III ENT model 40701620-1 is intended for use in the electromagnetic environment specified below.
 The customer or user of the UNIDRIVE® S III ENT model 40701620-1 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	0 V, 10 msec 92 V, 100 msec 168 V, 500 msec 0 V, 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Table 204
Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting

UNIDRIVE® S III ENT model **40701620-1** is intended for use in the electromagnetic environment specified below.
 The customer or user of the UNIDRIVE® S III ENT model **40701620-1** should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbance acc. to IEC 61000-4-6 Radiated RF disturbance acc. to IEC 61000-4-3	3 Vrms 150 kHz to < 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the UNIDRIVE® S III ENT model 40701620-1/40701720-1/40701820-1 , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3,5/3] \sqrt{P}$ 150 kHz to < 80 MHz $d = [3,5/3] \sqrt{P}$ 80 MHz to < 800 MHz $d = [3,5/3] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 206
Recommended separation distances between portable and mobile RF communications equipment
and the UNIDRIVE® S III ENT model 40701620-1

UNIDRIVE® S III ENT model 40701620-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot bear the cost of transportation or the risks associated with shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:
Company stamp/signature:

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KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen
Germany (Германия)

To be filled out by device owner:
Return address/company stamp:

Field of application:

Type of device: Serial no.:

Purchase date:

Signature/Date:

Please attach
sufficient
postage

Reply card

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen
Germany

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Absender/Firmenstempel:

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Geräte-Typ:

Serien-Nr.:

Kauf-Dat.:

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Antwortkarte

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen

KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postfach 230, 78503 Tuttlingen, Germany
Phone: +49 (0)7461 708-0,
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com

KARL STORZ Endoskope Berlin GmbH
Ohlauer Straße 43
10999 Berlin, Germany
Phone: +49 (0)30 30 69090
Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.
2345 Argenta Road, Suite 100
Mississauga, ON, L5N 8K4, Canada
Phone: +1 905 816-8100, Fax: +1 905 858-0933
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, +1 800 421-0837
Fax: +1 424 218-8526
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseas.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, Mexico
Phone: +52 5552 5056 07, Fax: +52 55 5545 0174

KARL STORZ Marketing América Do Sul Ltda.
R. Joaquim Floriano, 413 - 20º andar -
Itaim Bibi, São Paulo - SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Cervilio 4449 Piso 10º
1425 Buenos Aires C. F., Argentina
Phone: +54 11 4772 4545, Fax: +54 11 4772 4433
E-Mail: info@karlstorz.com.ar

KARL STORZ ENDOSKOPI NORGE AS
P.O. Box 153
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14,
12739 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Itälähdentkatu 23a, 00210 Helsinki, Finland
Phone: 0035-896824774, Fax: 0035-8968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopia France S. A.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Francia - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 071, Fax: +385 1 6406 077
E-Mail: b.vlajo@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: plsama@karlstorz.si

KARL STORZ Endoskope Greece M.E.P.E.
Sokratous & Kyprou 2, 151 27 Melissia, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

KARL STORZ Endoscopia ROMANIA srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

TOV KARL STORZ Ukraine
Prospekt Geroev Stalingrada, 18b
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19, -20
Fax: +380 44 42668-14

OOO KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Phoenicia Street, Mina Elhosn
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax: +961 1 365151

KARL STORZ Endoscopy South Africa (Pty) Ltd.
P.O. 6061,
Roggebaai
Cape Town 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSKOPY Kasachstan
Korgalynskoe shosse, 2A
Hotel „Duman“
Astana 010000, Kazakhstan
Phone/Fax: +7 7172 79 15 62, 79 15 63
E-Mail: info@karlstorz.kz

KARL STORZ Endoscopy - Gulf & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O.Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.ksegne@karlstorz-rce.com

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +84 991 8442, Fax: +84 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 - 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hol Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 (0)24 6258 9911
Fax: +86 (0)24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
#06-07 e-Centre @ Redhill
Singapore 159471
Phone: +65 65325548, Fax: +65 65323832
E-Mail: infoosp@karlstorz.com.sg
E-Mail: servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 (0)3 6380-8622,
Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com