Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Kathrine Mary Aaberg
December 5, 2019

Ashley Luloff
Sr Regulatory Affairs Specialist

December 5, 2019

Addition to Instructions for Use

1. Name of medical device	Kits of Leads for stimulation of the spinal cord and peripheral nerves in variants with accessories In the text may be indicated as "Lead", "Medical device", "Device", "MD".
2.1 Prescription of medical device	Spinal cord or peripheral nerve stimulation
2.2 Operational principles of medical device	The Medtronic SCS/PNS therapy systems are fully implantable and multiprogrammable systems. All of the Medtronic SCS/PNS systems are comprised of an implantable neurostimulator, a thin, insulated lead, and an extension to connect the neurostimulator and the lead. An external physician programmer is used to non-invasively adjust patient stimulation parameters. Patients may also make stimulation adjustments, based on their pain needs, via the external patient programmer. The neurostimulator portfolio includes 2 main families. The Restore family provides both non-rechargeable and fully rechargeable neurostimulators to accommodate high energy users. In addition, Itrel 4 is still offered as an option for low energy users. The Medtronic SCS/PNS therapy system provides a portfolio of leads for both percutaneous and surgical placement. Percutaneous leads are available with 4 or 8 electrode options while surgical leads have a 16 electrode configurations.
3. Information on potential consumers of medical device, use conditions	To be used in medical and preventive treatment facility. The lead should be used only by or under the supervision of physicians trained. Patient selection Select patients carefully to ensure that: ▪ their symptoms are of physiological origin. ▪ they are appropriate candidates for surgery.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	• they can properly operate the system. • they received satisfactory results from test stimulation.
4. Patient Counseling Information	Physicians should provide patients and caregivers with information about the following: • The components of the neurostimulation system: lead, extension, and neurostimulator • Instructions for using the neurostimulation system • The indications, contraindications, warnings, precautions, and adverse events for a neurostimulation system Physicians should also instruct patients to do the following: • Always inform any health care personnel that they have an implanted neurostimulation system before any procedure is begun. • Contact their physician if they notice any unusual symptoms or signs (for example, discomfort, paresthesia).
5. Description of the main functional parts of medical device	The Medtronic Models 3487A PISCES Quad, 3887 PISCES Quad Compact, or 3888 PISCES Quad Plus Leads The lead has electrodes on the distal end; the proximal end fits into a four conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead and performing intraoperative test stimulation.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com


Figure 1. Electrode, lead contact, and stylet. Model 3888 shown.

The Medtronic Models 3487A PISCES Quad, 3887 PISCES Quad Compact, or 3888 PISCES Quad Plus Leads

Description	Model 3487A	Model 3887	Model 3888
Connector	Quadrupolar in-line	Quadrupolar in-line	Quadrupolar in-line
Shape	Straight	Straight	Straight
Conductor resistance	<150 Ω (33 cm) (3487A-33) <200 Ω (45 cm) (3487A-45)	<150 Ω (33 cm) (3887-33) <200 Ω (45 cm) (3887-45)	<150 Ω (33 cm) (3888-33) <200 Ω (45 cm) (3888-45)

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	<250 Ω (56 cm) (3487A-56)	45)<250 Ω (56 cm) (3887-56)	<250 Ω (56 cm) (3888-56)
Length	20 cm – 110 cm	20 cm – 110 cm	20 cm – 110 cm
Diameter	1.3mm	1.3mm	1.3mm
Mass	1,04 g (33 cm) (3487A-33) 1,42 g (45 cm) (3487A-45) 1,8 g (56 cm) (3487A-56)	1,04 g (33 cm) (3887-33) 1,42 g (45 cm) (3887-45) 1,8 g (56 cm) (3887-56)	1,04 g (33 cm) (3888-33) 1,42 g (45 cm) (3888-45) 1,8 g (56 cm) (3888-56)
Distal end			
Number of electrodes	4	4	4
Electrode shape	Cylindrical	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm	3.0 mm
Electrode spacing	6.0mm	6.0mm	6.0mm
Electrode distance	30.0mm	24.0mm	60.0mm
Distal tip distance	Not applicable	Not applicable	7.9 mm
Radiopaque marker	No	No	Yes
Proximal end			
Lead contact length	2.3mm	2.3mm	2.3mm
Lead contact distance	15.2mm	15.2mm	15.2mm
Stylet handle length	40.0mm	40.0mm	40.0mm

The tolerances are $\pm 5\%$ if another value is not indicated

The Medtronic Vectris SureScan MRI 1x8 Subcompact Models 977A160, 977A175, and 977A190 Leads and the Vectris SureScan MRI 1x8 Compact Models 977A260, 977A275, and 977A290 Leads are part of a Neurostimulation systems with SureScan MRI Technology for pain therapy. The lead has electrodes on the distal end; the proximal (connector) end fits into an 8-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com


Figure 1. Electrode, lead contact, and stylet (shown with bent stylet).

Device specifications for leads: Models 977A1 and 977A2 series

Description	Models 977A160, 977A175, 977A190	Models 977A260, 977A275, 977A290
Connector	Octapolar, in-line	Octapolar, in-line
Conductor resistance	1.6 Ω per cm	1.6 Ω per cm
Length	60, 75, 90 cm	60, 75, 90 cm
Diameter(lead body)	1.3mm	1.3mm

Medtronic

Medtronic


710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Mass	2,8 g (977A160) 3,5 g (977A175) 4,2 g (977A190)	2,8 g (977A260) 3,5 g (977A275) 4,2 g (977A290)
Distal(electrode) end		
Number of electrodes	8	8
Electrode shape	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm
Electrode spacing	1.5mm	1.5mm
Electrode distance	34.5mm	52.0mm
Distal tip distance	Not applicable	Not applicable
Radiopaque tip marker	No	No
Proximal (connector) end		
Lead contact length	1.5mm	1.5mm
'O'contact length	3.0mm	3.0mm
Lead contact spacing	2.8mm	2.8mm
Lead contact distance	22.5mm	22.5mm
Stylet handle length	29.0mm	29.0mm
Stylet wire diameter	0.25, 0.30 mm	0.25, 0.30 mm
The tolerances are $\pm 5\%$ if another value is not indicated		
The lead Vectris 1x8 compact 977D160, Vectris 1x8 subcompact 977D260 has electrodes on the distal end; the proximal (connector) end fits into an 8-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead.		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Distal **Proximal**
Electrode, lead contact, and stylet (shown with bent stylet).

Device specifications for leads: Models 977D160 and 977D260

Description	Models 977D160	Models 977D260
Connector	Octapolar, in-line	Octapolar, in-line
Conductor resistance	1.6 Ω per cm	1.6 Ω per cm
Length	60cm	60cm
Diameter (lead body)	1.3 mm	1.3 mm
Mass	2,8g	2,8g
Distal (electrode) end		
Number of electrodes	8	8
Electrode shape	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Electrode spacing	1.5mm	4.0mm
Electrode distance	34.5 mm	52.0 mm
Distal tip distance	Not applicable	Not applicable
Radiopaque tip marker	no	no
Proximal (connector) end		
Lead contact length	1.5 mm	1.5 mm
'0' contact length	3.0 mm	3.0 mm
Lead contact spacing	2.8 mm	2.8 mm
Lead contact distance	22.5 mm	22.5 mm
Stylet handle length	29.0 mm	29.0 mm
Stylet wire diameter	0.25 mm, 0.30 mm	0.25 mm, 0.30 mm

The tolerances are $\pm 5\%$ if another value is not indicated


The Medtronic Specify SureScan MRI Surgical Leads, Models 977C165, 977C190, 977C265, and 977C290, are part of a neurostimulation system for pain therapy.

The lead has 16 electrodes on the distal end; the two proximal ends fit into 8-conductor connectors (ie, a Medtronic neurostimulator). One lead body has a white marker band to indicate electrodes 0-7.


Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Specify SureScan MRI 5-6-5 Models 977C165 and 977C190 Surgical Lead (electrode side shown).

Specify SureScan MRI 2x8 Models 977C265 and 977C290 Surgical Lead (electrode side shown).

Device specifications:

Description	Models 977C165, 977C190	Models 977C265, 977C290
Connector	Octapolar, in-line	Octapolar, in-line
Shape	Contoured	Contoured
Conductor resistance	<77 Ω (65 cm) <107 Ω (90 cm)	<77 Ω (65 cm) <107 Ω (90 cm)
Length (includes paddle)	65 cm, 90 cm	65 cm, 90 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Diameter	1.3 mm	1.3 mm
	Mass	4,39 g (65 cm)	4,39 g (65 cm)
		6,08 g (90 cm)	6,08 g (90 cm)
	Distal (electrode) end		
	Number of electrodes	16	16
	Electrode shape	Rectangular	Rectangular
	Electrode size (width x length)	1.5 mm x 4.0 mm	1.5 mm x 4.0 mm
	Electrode stimulating area	6.0 mm ²	6.0 mm ²
	Electrode spacing (edge to edge)		
	in-line spacing	4.5 mm	1.0 mm
	row spacing	1.0 mm	1.0 mm
	Electrode paddle length	64.2 mm	56.4 mm
	Electrode paddle width	10 mm	7.6 mm
	Electrode paddle thickness	2.0 mm	1.9 mm
	Proximal (connector) end		
	Lead contact length	1.5 mm	1.5 mm
	Lead contact spacing (center to center)	2.8 mm	2.8 mm
	Lead contact distance	22.5 mm	22.5 mm
	The tolerances are $\pm 5\%$ if another value is not indicated		
	Accessories of leads		

Quadripolar extension (8-4), model 37083-60, 37083-40, 37083-20
The extension has four connectors on the distal end and eight connectors on the proximal end. The in-line distal end connects to a Medtronic lead. The proximal end connects to a Medtronic neurostimulator.

Medtronic


Medtronic

710 Medtronic Parkway

Minneapolis MN 55432-5604

USA

www.medtronic.com

Device specifications for the Model 37083 Extension

Description	Model 37083
Resistance	Maximum 24.0 Ω for all lengths
Length	10 cm – 110 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Mass	4,32 g (60 cm) 2,88 g (40 cm) 1,44 g (20 cm)
Distal end (lead)	
Connector	Quadripolar, in-line
Contactspacing	4.3 mm
Leadentrancediameter	1.5 mm
Outerdiameter	3.8 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Contactspacing	2.8 mm
Diameter	1.3 mm

The tolerances are $\pm 5\%$ if another value is not indicated

Dual quadripolar extension (8-2-4), model 37082-60, 37082-40, 37082-20

The extension has eight connectors on the distal end, four connectors on each distal leg, and eight connectors on the proximal end. The in-line distal end is bifurcated to connect to one or two Medtronic quadripolar leads. The proximal end connects to a Medtronic neurostimulator.

Note: The marker band is for lead identification. The lead connected to the extension leg with the marker band is for

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

either circuits 0-3 or circuits 8-11.

Device specifications for the Model 37082 Extension

Description	Model 37082
Resistance	Maximum 24.0 Ω for all lengths
Length	10 cm – 110 cm
Mass	9,32 g (60 cm) 6,21 g (40 cm) 3,1 g (20 cm)
Distal end (lead)	
Connector	Quadripolar, in-line
Contact spacing	4.3 mm
Lead entrance diameter	1.5 mm
Outer diameter	3.8 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Contact spacing	2.8 mm
Diameter	1.3 mm

The tolerances are $\pm 5\%$ if another value is not indicated


1x8 Extension, model 37081-60, 37081-40, 37081-20

The extension has eight connectors each on the distal end and the proximal end. The in-line distal end connects to a Medtronic lead. The proximal end connects to a Medtronic neurostimulator.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Device specifications for the Model 37081 Extension

Description	Model 37081
Resistance	Maximum 16.0 Ω for all lengths
Length	10 cm – 110 cm
Contact spacing	2.8mm
Mass	4,26 g (60 cm) 2,84 g (40 cm) 1,42 g (20 cm)
Distal end (lead)	
Connector	Quadripolar, in-line
Conductor wire diameter	1.3 mm
Lead entrance diameter	4.7 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Diameter	1.3 mm

The tolerances are $\pm 5\%$ if another value is not indicated

Multi-lead Trialing Cable, model 355531


The multi-lead trialing cable has a lead housing with a lockable lid on the distal end and a connector plug on the proximal end (Figure 1). The lead housing has channels for up to four quadripolar (4-contact) leads and two octapolar (8-contact)

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

leads. The connector plug inserts into an output jack on an external neurostimulator (ENS).

Device specifications for the Multi-lead Trialing Cable Model 355531

Description	Value
Leads	4 quadripolar; 2 octapolar
Electrodes	32 electrodes, supporting 16 active
Mass	50 g
Width (housing and slide lock)	45 mm
Length (housing)	57 mm
Thickness (housing)	10 mm
Volume (housing)	18 cm ³
Lid	Clear, with lid lift tab
Latch type	Snap lid with slide lock
Length (housing, cable, and plug)	152 cm
Storage temperature	-34° to +57°C

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com


Models 3487A PISCES Quad, 3887, 3888 PISCES Quad Compact PISCES Quad Plus Leads

Percutaneous extension

Full length	13+/-1 in / 33,02+/-2,54 cm
Tubing tip length	0,826+/-0,2 in / 2,10+/-0,05 cm
Mass	1,77 g

The tolerances are $\pm 5\%$ if another value is not indicated

Drawing of percutaneous extension


Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA


www.medtronic.com

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

EZ Anchor

Description	Value
Inner diameter	0,12 in / 0,3 cm
Outer diameter	0,14 in / 0,36 cm
Total length	2,25+/-0,1 in / 5,71+/-0,25 cm
Mass	0,5 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

2 Wing Anchor

Full length	0,6+/-0,02 in / 1,52+/-0,05 cm
Outer diameter	0,12 in / 0,33 cm
Inner diameter	0,02 in / 0,05 cm
Mass	0,15 g

The tolerances are $\pm 5\%$ if another value is not indicated

3 Wing Anchor

Full length	0,81+/-0,02 in / 2,06+/-0,05 cm
Outer diameter	0,13 in / 0,33 cm
Inner diameter	0,058+/-0,03 in / 0,15+/-0,04 cm
Mass	0,2 g

The tolerances are $\pm 5\%$ if another value is not indicated

Twist Lock anchor

Dimensions	22,35 mm x 10,21 mm x 10 mm
Mass	0,72 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Connector boot

Outer diameter	0,195+/-0,003 in / 0,495+/-0,007 cm
Outer tip diameter	0,024 in / 0,061 cm
Base length	1,096+0,015+0,01 in / 2,784+0,038-0,025 cm
Tip length	0,22+/-0,005 in / 0,559+/-0,127 cm
Mass	0,38 g

The tolerances are $\pm 5\%$ if another value is not indicated

Guide wire

Length	28 cm
Outer diameter	0,84 mm
Mass	0,76 g

The tolerances are $\pm 5\%$ if another value is not indicated

Needle, 15 gauge (epidural, curved or modified Tuohy)

Medtronic

Medtronic


710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Assy length	110,95 mm
	Full length of needle	105,48 mm
	Working length of needle	86,45 mm
	Outer diameter of needle	1,750 mm - 1,900 mm
	Minimum inner diameter of needle	1,3 mm
	Working length of mandrin	93,75 mm
	Outer diameter of mandrin	0,87 mm
	Hardness	≤92 HRC
	Roughness	Ra≤0,16 μm
	Bending strength of needle	25 N ± 0,1 N
	Blunting radius	≤ 0,03 mm
	Mass	5,9 g
	Luer adaptor	Luer-lock
	The tolerances are ±5% if another value is not indicated	
Drawing of epidural needle		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Screener cable

Cable length	146 cm
Cable outer diameter	2,96 cm
Mass	21,04 g

The tolerances are $\pm 5\%$ if another value is not indicated

Straight stylet 0,25 mm

Wire length	30,19 $\pm$ 0,51 cm
Mass	0,41 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA


www.medtronic.com

Short stylet 0,25 mm

Wire length	2,58 mm
Mass	0,23 g

The tolerances are $\pm 5\%$ if another value is not indicated

Drawing of the stylet (for straight stylet as example)

Tunneling Rod

Tubing diameter	3,2 mm
-----------------	--------

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Outer diameter of mounting thread	3,11 mm
Tubing length	8,14 in / 20,68 cm
Mass	21,26 g
Hardness	≤92 HRC
Roughness	Ra≤3,2 μm

The tolerances are ±5% if another value is not indicated

Tunneling Tube

Outer diameter	0,192+/-0,005 in / 4,9+/-0,1 mm
Inner diameter	2,62 mm

The tolerances are ±5% if another value is not indicated

Tunneling Tip

Diameter	3,92 mm
Inner diameter of mounting thread	2,6 mm
Tip length	1,01+/-0,06 in / 2,56 +/-0,15 cm
Tip angle	15+/-3 °
Blunting radius	≤0,03 mm
Mass (with tunneling tube)	1,7 g
Hardness	≤92 HRC
Roughness	Ra≤3,2 μm

The tolerances are ±5% if another value is not indicated

Medtronic


Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA

www.medtronic.com

Drawings of tunneling instruments: tunneling rod, tunneling tip, tunneling tube


Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA

www.medtronic.com



Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA

www.medtronic.com

Wrench, torque

Medtronic


Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Hex tip outer diameter	From 0,40 in to 0,45 in / 0,10 to 0,11 cm
Screw outer diameter	0,434 in / 1,1 cm
Screw length	1,15+/-0,15 in / 2,92+/-0,38 cm
Hex tip length	,055+/-0,15 in / 1,40+/-0,38 cm
Mass	2,93 g
Hardness	48 HRC
Roughness	Ra≤3,2 μm

The tolerances are ±5% if another value is not indicated

Drawing of torque wrench

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Models 977A160, 977A175, and 977A190 Leads and the Vectris SureScan MRI 1x8 Compact Models 977A260, 977A275, and 977A290 Leads

Anchor with anchor dispenser tool

Full height	128,83 mm
Tip height	1,32 in / 3,35 cm
Tip outer diameter	0,188 in / 0,48 cm
Shaft outer diameter	0,35 in / 0,89 cm
Handle outer diameter	0,45 in / 1,09 cm
Mass of assy	10,13 g

The tolerances are $\pm 5\%$ if another value is not indicated

Spare anchor on spare handle

Full height	2,42 in / 6,15 cm
Tip height	1,32 in / 3,35 cm
Tip outer diameter	0,188 in / 0,48 cm
Shaft outer diameter	0,35 in / 0,89 cm
Handle outer diameter	0,45 in / 1,09 cm
Mass of assy	4,5 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Anchor removal tool

Full length	1,62+/-0,05 in / 4,11+/-0,13 cm
Full height	0,71+/-0,05 in / 1,80+/-0,13 cm
Mass	1,33 g

The tolerances are $\pm 5\%$ if another value is not indicated

Guide wire

Length	28 cm
Outer diameter	0,84 mm
Mass	0,76 g

The tolerances are $\pm 5\%$ if another value is not indicated

Needle – curved tip, 14 gauge, 8.9 cm

Assy length	114,25 mm
Total length of needle	108,78 mm
Working length of needle	90,17 mm
Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,6 mm
Working length of mandrin	102,56 mm
Outer diameter of mandrin	1,39 mm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Hardness	≤92 HRC
Roughness	Ra≤0,16 μm
Bending strength of needle	40 N ± 0,1 N
Blunting radius	≤ 0,03 mm
Mass	6,48 g
Luer adaptor	Luer-lok

The tolerances are ±5% if another value is not indicated

Needle – straight tip, 14 gauge, 8.9 cm

Assy length	111,95 mm
Total length of needle	105,97 mm
Working length of needle	87,36 mm
Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,6 mm
Working length of mandrin	99,44 mm
Outer diameter of mandrin	1,23 mm
Hardness	≤92 HRC
Roughness	Ra≤0,16 μm
Bending strength of needle	40 N ± 0,1 N
Blunting radius	≤ 0,03 mm
Mass	6,38 g

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Luer adaptor	Luer-lok
The tolerances are $\pm 5\%$ if another value is not indicated	
Stylet – straight, 0.30 mm diameter (green stylet handle and green cap)	
Wire length	30,19 $\pm$ 0,51 cm
Mass	0,89 g
The tolerances are $\pm 5\%$ if another value is not indicated	
Stylet – bent, 0.25 mm diameter (white stylet handle and blue cap)	
Wire length	30,19 $\pm$ 0,51 cm
Mass	0,82 g
The tolerances are $\pm 5\%$ if another value is not indicated	
Stylet – straight, 0.25 mm diameter (white stylet handle and white cap)	
Wire length	30,19 $\pm$ 0,51 cm
Mass	0,8 g
The tolerances are $\pm 5\%$ if another value is not indicated	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Tunneling Rod

Tubing diameter	3,2 mm
Outer diameter of mounting thread	3,11 mm
Tubing length	8,14 in / 20,68 cm
Mass	21,26 g
Hardness	≤92 HRC
Roughness	Ra≤3,2 μm

The tolerances are ±5% if another value is not indicated

Tunneling Tube

Outer diameter	0,192+/-0,005 in / 4,9+/-0,1 mm
Inner diameter	2,62 mm

The tolerances are ±5% if another value is not indicated

Tunneling Tip

Diameter	3,92 mm
Inner diameter of mounting thread	2,6 mm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Tip length	1,01+/-0,06 in / 2,56 +/-0,15 cm
Tip angle	15+/-3 °
Blunting radius	≤0,03 mm
Mass (with tunneling tube)	1,7 g
Hardness	≤92 HRC
Roughness	Ra≤3,2 μm

The tolerances are ±5% if another value is not indicated

Wrench, torque

Hex tip outer diameter	From 0,40 in to 0,45 in / 0,10 to 0,11 cm
Screw outer diameter	0,434 in / 1,1 cm
Screw length	1,15+/-0,15 in / 2,92+/-0,38 cm
Hex tip length	,055+/-0,15 in / 1,40+/-0,38 cm
Mass	2,93 g
Hardness	48 HRC
Roughness	Ra≤3,2 μm

The tolerances are ±5% if another value is not indicated

Models 977D160 and the Vectris 1x8 Compact Models 977D260 Leads

Guide wire

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Length	28 cm
Outer diameter	0,84 mm
Mass	0,76 g

The tolerances are $\pm 5\%$ if another value is not indicated

Needle – curved tip, 14 gauge, 8.9 cm

Assy length	114,25 mm
Total length of needle	108,78 mm
Working length of needle	90,17 mm
Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,6 mm
Working length of mandrin	102,56 mm
Outer diameter of mandrin	1,39 mm
Hardness	≤ 92 HRC
Roughness	$Ra \leq 0,16 \mu m$
Bending strength of needle	$40 N \pm 0,1 N$
Blunting radius	$\leq 0,03 mm$
Mass	6,48 g
Luer adaptor	Luer-lok

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Needle – straight tip, 14 gauge, 8.9 cm

Assy length	111,95 mm
Total length of needle	105,97 mm
Working length of needle	87,36 mm
Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,6 mm
Working length of mandrin	99,44 mm
Outer diameter of mandrin	1,23 mm
Hardness	≤92 HRC
Roughness	Ra≤0,16 µm
Bending strength of needle	40 N ± 0,1 N
Blunting radius	≤ 0,03 mm
Mass	6,38 g
Luer adaptor	Luer-lok

The tolerances are ±5% if another value is not indicated

Stylet – straight, 0.30 mm diameter (green stylet handle and green cap)

Wire length	30,19+/-0,51 cm
Mass	0,89 g

The tolerances are ±5% if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Stylet – bent, 0.25 mm diameter (white stylet handle and blue cap)

Wire length	30,19+/-0,51 cm
Mass	0,82 g

The tolerances are $\pm 5\%$ if another value is not indicated

Stylet – straight, 0.25 mm diameter (white stylet handle and white cap)

Wire length	30,19+/-0,51 cm
Mass	0,8 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic Specify SureScan MRI Surgical Leads, Models 977C165, 977C190, 977C265, and 977C290 Package contents:

EZ Anchor

Description	Value
Inner diameter	0,12 in / 0,3 cm
Outer diameter	0,14 in / 0,36 cm
Total length	2,25+/-0,1 in / 5,71+/-0,25 cm
Mass	0,5 g

The tolerances are $\pm 5\%$ if another value is not indicated

2 Wing Anchor

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	0,6+/-0,02 in / 1,52+/-0,05 cm
Outer diameter	0,12 in / 0,33 cm
Inner diameter	0,02 in / 0,05 cm
Mass	0,15 g

The tolerances are $\pm 5\%$ if another value is not indicated

3 Wing Anchor

Full length	0,81+/-0,02 in / 2,06+/-0,05 cm
Outer diameter	0,13 in / 0,33 cm
Inner diameter	0,058+/-0,03 in / 0,15+/-0,04 cm
Mass	0,2 g

The tolerances are $\pm 5\%$ if another value is not indicated

Passing elevator

Tip height	2,06+/-0,1 in / 5,23+/-0,25 cm
Tip width	1,82+/-0,2 in / 4,62+/-0,51 cm
Tip outer diameter	0,3 in / 0,76 cm
Tip inner diameter	0,088+/-0,1 in / 0,22+/-0,02 cm
Mass	3,04 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Lead blank

Dimensions	11,5 cm x 0,7 cm
Mass	1,42 g

The tolerances are $\pm 5\%$ if another value is not indicated

Tunneling Rod

Tubing diameter	3,2 mm
Outer diameter of mounting thread	3,11 mm
Tubing length	8,14 in / 20,68 cm
Mass	21,26 g
Hardness	≤ 92 HRC
Roughness	$Ra \leq 3,2 \mu m$

The tolerances are $\pm 5\%$ if another value is not indicated

Tunneling Tube

Outer diameter	0,192 $\pm$ 0,005 in / 4,9 $\pm$ 0,1 mm
Inner diameter	2,62 mm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

The tolerances are $\pm 5\%$ if another value is not indicated

Tunneling Tip

Diameter	3,92 mm
Inner diameter of mounting thread	2,6 mm
Tip length	1,01 $\pm$ 0,06 in / 2,56 $\pm$ 0,15 cm
Tip angle	15 $\pm$ 3 °
Blunting radius	$\leq$ 0,03 mm
Mass (with tunneling tube)	1,7 g
Hardness	$\leq$ 92 HRC
Roughness	Ra $\leq$ 3,2 μ m

The tolerances are $\pm 5\%$ if another value is not indicated

Wrench, torque

Hex tip outer diameter	From 0,40 in to 0,45 in / 0,10 to 0,11 cm
Screw outer diameter	0,434 in / 1,1 cm
Screw length	1,15 $\pm$ 0,15 in / 2,92 $\pm$ 0,38 cm
Hex tip length	,055 $\pm$ 0,15 in / 1,40 $\pm$ 0,38 cm
Mass	2,93 g
Hardness	48 HRC
Roughness	Ra $\leq$ 3,2 μ m

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

The tolerances are $\pm 5\%$ if another value is not indicated

1x8 Extension, model 37081-60, 37081-40, 37081-20

Connector boot, 1 x 8 low-profile

Length	2,3 cm
Diameter	0,5 cm
Mass	0,18 g

The tolerances are $\pm 5\%$ if another value is not indicated

Set screw

Length	1,35 mm
Outer diameter of thread	1,77 mm
Inner diameter of thread	1,67 mm
Thread pitch	0,34 mm
Mass	0,0114 g

The tolerances are $\pm 5\%$ if another value is not indicated

Tunneling tools: extension passer and obturator

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Length	39 cm
	Handle full width	0,22 in / 0,56 cm
	Handle tip length	0,265+0,01-0,005 in / 0,67+0,02-0,01 cm
	Handle outer diameter	0,17+0,003-0,12 in / 0,43+0,07-0,03 cm
	Tubing outer diameter	0,086+/-0,006 in / 0,22+/-0,01 cm
	Tip smaller outer diameter	0,057+/-0,003 in / 0,145+/-0,007 cm
	Tip larger outer diameter	0,115+/-0,003 in / 0,292+/-0,018 cm
	Tip widening angle	30 °
	Mass with handle	31.88 g
	Hardness	≤92 HRC
	Roughness	Ra≤3,2 μm
	The tolerances are ±5% if another value is not indicated	
	Wrench, torque	
	Hex tip outer diameter	From 0,40 in to 0,45 in / 0,10 to 0,11 cm
	Screw outer diameter	0,434 in / 1,1 cm
	Screw length	1,15+/-0,15 in / 2,92+/-0,38 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Hex tip length	,055+/-0,15 in / 1,40+/-0,38 cm
	Mass	2,93 g
	Hardness	48 HRC
	Roughness	Ra≤3,2 μm
	The tolerances are ±5% if another value is not indicated	
	Models 37082, 37083 Extension	
	Connector boot, 1 x 4 low-profile	
	Length	3,4 cm
	Diameter	0,5 cm
	Mass	0,27 g
	The tolerances are ±5% if another value is not indicated	
	Octapolar in-line neurostimulator plug	
	Dimensions	2,4 cm x 0,7 cm
	Plug diameter	1 mm
	Mass	0,11 g
	The tolerances are ±5% if another value is not indicated	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Set screw

Length	1,35 mm
Outer diameter of thread	1,77 mm
Inner diameter of thread	1,67 mm
Thread pitch	0,34 mm
Mass	0,0114 g

The tolerances are $\pm 5\%$ if another value is not indicated

Tunneling tools: extension passer and obturator

Length	39 cm
Handle full width	0,22 in / 0,56 cm
Handle tip length	0,265 $\pm$ 0,01-0,005 in / 0,67 $\pm$ 0,02-0,01 cm
Handle outer diameter	0,17 $\pm$ 0,003-0,12 in / 0,43 $\pm$ 0,07-0,03 cm
Tubing outer diameter	0,086 $\pm$ -0,006 in / 0,22 $\pm$ -0,01 cm
Tip smaller outer diameter	0,057 $\pm$ -0,003 in / 0,145 $\pm$ -0,007 cm
Tip larger outer diameter	0,115 $\pm$ -0,003 in / 0,292 $\pm$ -0,018 cm
Tip widening angle	30 °
Mass with handle	31.88 g
Hardness	≤ 92 HRC

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Roughness		Ra≤3,2 μm
The tolerances are ±5% if another value is not indicated		
Carrier		
Full length		19,24 in / 48,87 cm
Tubing length		16,69 in / 42,39 cm
Handle outer diameter		0,23+/-0,003 in / 0,584+/-0,008 cm
Handle inner diameter		0,185 in / 0,47 cm
Mass		1,91 g
The tolerances are ±5% if another value is not indicated		
Wrench, torque		
Hex tip outer diameter		From 0,40 in to 0,45 in / 0,10 to 0,11 cm
Screw outer diameter		0,434 in / 1,1 cm
Screw length		1,15+/-0,15 in / 2,92+/-0,38 cm
Hex tip length		,055+/-0,15 in / 1,40+/-0,38 cm
Mass		2,93 g
Hardness		48 HRC
Roughness		Ra≤3,2 μm
The tolerances are ±5% if another value is not indicated		


Boots and anchors accessory kit, model 3550-24

The Medtronic Model 3550-24 Accessory Kit includes boots and anchors that are part of a neurostimulation system for pain therapy.

Anchors provide a clinician with an anchoring solution with a holding strength between the anchor and lead that is not dependent on suturing technique or the use of medical adhesive.

Boots are used to cover the distal end of an extension. The purpose of the boot is to prevent the ingress of body fluids into the extension in the event only one bore is used.

EZ Anchor

Description	Value
Inner diameter	0,12 in / 0,3 cm

Medtronic

Medtronic


710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Outer diameter	0,14 in / 0,36 cm
	Total length	2,25+/-0,1 in / 5,71+/-0,25 cm
	Mass	0,5 g
	The tolerances are ±5% if another value is not indicated	
2 Wing Anchor		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	0,6+/-0,02 in / 1,52+/-0,05 cm
Outer diameter	0,12 in / 0,33 cm
Inner diameter	0,02 in / 0,05 cm

Medtronic

Medtronic


710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Mass

0,15 g

The tolerances are $\pm 5\%$ if another value is not indicated

3 Wing Anchor

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	0,81+/-0,02 in / 2,06+/-0,05 cm
Outer diameter	0,13 in / 0,33 cm
Inner diameter	0,058+/-0,03 in / 0,15+/-0,04 cm
Mass	0,2 g

The tolerances are $\pm 5\%$ if another value is not indicated

Low Profile 1x4 Boot

Outer diameter	0,195+/-0,003 in / 0,495+/-0,007 cm
Outer tip diameter	0,024 in / 0,061 cm
Base length	1,096+0,015+0,01 in / 2,784(+0,038/-0,025) cm
Tip length	0,22+/-0,005 in / 0,559+/-0,127 cm
Mass	white: 0,27 g transparent: 0,24 g

The tolerances are $\pm 5\%$ if another value is not indicated

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Low Profile 1 x 8 Closed Boot



Full length	0,688+/-0,02 in / 1,75+/-0,05 cm
Inner diameter	0,168 in / 0,43 cm
Outer diameter	0,21 in / 0,53 cm
Width	0,23 in / 0,58 cm
Mass	0,15 g

The tolerances are $\pm 5\%$ if another value is not indicated

Flat Boot

Medtronic

Medtronic
710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	0,7 in / 1,78 cm
Inner height	0,52 in / 1,32 cm
Outer width	0,18 in / 0,46 cm
Inner width	0,13 in / 0,33 cm
Mass	0,42 g


Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

The tolerances are $\pm 5\%$ if another value is not indicated

1x8 Boot

Full length	1,19+0,01-0,02 in / 3,02(+0,02/-0,05) cm
Outer diameter	6,2 mm
Inner diameter	5,9 mm
Base length	2,4 cm
Mass	White: 0,4 g Transparent: 0,37 g


The tolerances are $\pm 5\%$ if another value is not indicated

1x4 Boot

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	2,3 cm
Outer diameter	4,8 mm
Inner diameter	4,5 mm
Base length	1,7 cm
Mass	White: 0,2 g Transparent: 0,18g

The tolerances are $\pm 5\%$ if another value is not indicated

Low profile 1x8 Boot

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	4,8 cm
Outer diameter	6,5 mm
Inner diameter	5,85 mm
Base length	4,15 cm
Mass	0,61 g

The tolerances are $\pm 5\%$ if another value is not indicated

Injex , models 97791, 97792

The Injex Model 97791(Bumpy), 97792(Bi-wing) Anchor Accessory Kit is designed only for use with Medtronic percutaneous in-line connector neurostimulation leads for spinal cord stimulation (SCS). The tools in this accessory kit are designed solely for the placement and removal of the Injex anchors.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Anchor dispenser tool

Full height	128,83 mm
Tip height	1,32 in / 3,35 cm
Tip outer diameter	0,188 in / 0,48 cm
Shaft outer diameter	0,35 in / 0,89 cm
Handle outer diameter	0,45 in / 1,09 cm
Mass of assy	10,13 g

The tolerances are $\pm 5\%$ if another value is not indicated

Anchor Removal Tool (ART)

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	1,62+/-0,05 in / 4,11+/-0,13 cm
Full height	0,71+/-0,05 in / 1,80+/-0,13 cm
Mass	1,33 g

The tolerances are $\pm 5\%$ if another value is not indicated

Bumpy Anchor 97791


Length	$2,71018 \pm 0,0762$ cm
Inner diameter	$0,08128 \pm 0,0127$ cm
Outer diameter	$0,3175 \pm 0,0127$ cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Bi-wing Anchor 97792	
	Length	2,71018 ± 0,0762 cm
	Inner diameter	0,08128 ± 0,0127 cm
	Outer diameter	0,3175 ± 0,0127 cm
 Inserting lead into ADT		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA

www.medtronic.com


Figure 8. Anchor Removal Tool


Model 3550-38 1x 8 Tunneling Tool Accessory Kit

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Handle, Obturator

Width	2,46 in / 6,24 cm
Diameter	0,228 in / 0,58 cm
Handle height	0,79 in / 2,01 cm

The tolerances are $\pm 5\%$ if another value is not indicated


Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA

www.medtronic.com

DETAIL A

DETAIL B

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Length	44,4 cm
	Handle full width	0,22 in / 0,56 cm
	Handle tip length	0,265+0,01-0,005 in / 0,67(+0,02/-0,01)cm
	Handle outer diameter	0,17+0,003-0,12 in / 0,43(+0,07/-0,03) cm
	Tubing outer diameter	0,086+/-0,006 in / 0,22+/-0,01 cm
	Tip smaller outer diameter	0,057+/-0,003 in / 0,145+/-0,007 cm
	Tip larger outer diameter	0,115+/-0,003 in / 0,292+/-0,018 cm
	Tip widening angle	30 °
	Mass with handle	32,1 g
	Hardness	≤92 HRC
	Roughness	Ra≤3,2 μm
	The tolerances are ±5% if another value is not indicated	
Carrier		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length	19,24 in / 48,87 cm
Tubing length	16,69 in / 42,39 cm
Handle outer diameter	0,23+/-0,003 in / 0,584+/-0,008 cm
Handle inner diameter	0,185 in / 0,47 cm
Mass	2,1 g

The tolerances are $\pm 5\%$ if another value is not indicated

Low profile 1x8 boot

Full length	4,8 cm
Outer diameter	6,5 mm
Inner diameter	5,85 mm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Base length	4,15 cm
	Mass	0,61 g
	The tolerances are ±5% if another value is not indicated	
Model 3550-43 Epidural Needle Accessory Kit		


Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA

www.medtronic.com

Full length	6,96-0,12 in / 17,68-0,30 cm
Working length of needle	15,0 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,5 mm
Working length of mandrin	16,4 cm
Outer diameter of mandrin	1,42 mm
Hardness	≤92 HRC
Roughness	Ra≤0,16 μm
Bending strength of needle	40 N ± 0,1 N
Blunting radius	≤ 0,03 mm
Mass	7,76 g
Luer adaptor	Luer-lok

The tolerances are ±5% if another value is not indicated

Model 3550-60 1x 8, Model 3555-60 Tunneling Tool Accessory Kit extension passer

The Medtronic® Model 3555-60 Extension Passer Accessory Kit is designed for use with the Medtronic Neurological Extensions. The accessory kit contains one extension passer, one obturator, and one carrier. These are sterile, single-use devices.

The Model 3555-60 Extension Passer Accessory Kit is indicated for Spinal Cord Stimulation (SCS) and Peripheral Nerve

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Stimulation (PNS) to aid in the treatment of chronic intractable pain of the trunk or limbs.

Product Specifications:

- The extension passer handle will require a minimum of 10 pounds of force to separate from the stainless steel tubing.
- The stainless steel tubing shall be designed to withstand one 90 degree bend without kinking. A Model 7489 extension shall be able to pass through the tube after being bent to 90 degrees. The extension passer may be straightened to an angle greater than 90 degrees to facilitate passage.
- The carrier tip shall be designed to resist breakage. When installed in the extension passer, the tip shall be able to be bent once to a 45 degree angle without permanent functional damage.
- The Extension Passer System shall meet the specification requirements over the operating temperature range of 68°F to 105°F.
- The Extension Passer System shall show no signs of damage or degradation by exposure to the process temperature extremes of -40°F to 167°F.
- The extension Passer System, as packaged, must withstand a sinusoidal vibration sweep from 5 Hz to 500 Hz to 5 Hz in each of three orthogonal axes and still meet the requirements of this specification.
- The Extension Passer System, as packaged, shall meet the requirements of this specification after being subjected to the mechanical shock associated with a 48-inch drop to a hard surface.
- The packaging and contents shall show no sign of damage nor degradation by exposure to storage temperature extremes of -30°F to 135°F (-34°C to 57°C).

Handle, Obturator

Width	2,46 in / 6,24 cm
Diameter	0,228 in / 0,58 cm
Handle height	0,79 in / 2,01 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

The tolerances are $\pm 5\%$ if another value is not indicated

Length	25,21 in / 64,03 cm
Handle full width	0,22 in / 0,56 cm
Handle tip length	0,265+0,01-0,005 in / 0,67(+0,02/-0,01) cm
Handle outer diameter	0,17+0,003-0,012 in / 0,43(+0,07/-0,03) cm
Tubing outer diameter	0,086+/-0,006 in / 0,22+/-0,01 cm
Tip smaller outer diameter	0,057+/-0,003 in / 0,145+/-0,007 cm
Tip larger outer diameter	0,115+/-0,007 in / 0,292+/-0,018 cm
Tip widening angle	30 °
Mass with handle	53 g
Hardness	≤ 92 HRC
Roughness	Ra $\leq 3,2$ μm

The tolerances are $\pm 5\%$ if another value is not indicated

Carrier

Medtronic

Medtronic


710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Full length	27,88 in / 70,81 cm
	Tubing length	25,33 in / 64,34 cm
	Handle outer diameter	0,23+/-0,003 in / 0,584+/-0,008 cm
	Handle inner diameter	0,185 in / 0,47 cm
	Mass	2,79 g
The tolerances are $\pm 5\%$ if another value is not indicated		
Low profile 1x8 boot		
	Full length	4,8 cm
	Outer diameter	6,5 mm
	Inner diameter	5,85 mm
	Base length	4,15 cm
	Mass	0,61 g
The tolerances are $\pm 5\%$ if another value is not indicated		
Model 3550-P4 Passing Elevator Accessory Kit		
Passing Elevator 1		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604

USA


www.medtronic.com

	Tip width	1,82+/-0,2 in / 4,75+/-0,51 cm
	Mass	2,25 g
	The tolerances are ±5% if another value is not indicated	
Passing Elevator 2		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Full length

8,56+/-0,25 in / 21,08+/-0,63 cm

Medtronic

Medtronic


710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Tip outer diameter	0,3 in / 0,76 cm
	Tip length	2,75+/-0,06 in / 6,86+/-0,15 cm
	Tubing outer diameter	0,125 in / 0,31 cm
	Elevating tip height	1,62+/-0,03 in / 4,11+/-0,08 cm
	Elevating tip width	1,5+/-0,09 in / 3,81+/-0,23 cm
	Mass	2,32 g
The tolerances are $\pm 5\%$ if another value is not indicated		
Passing Elevator 3		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com


SECTION A-A

Tip height	2,06+/-0,1 in / 5,23+/-0,25 cm
Tip width	1,82+/-0,2 in / 4,62+/-0,51 cm
Tip outer diameter	0,3 in / 0,76 cm
Tip inner diameter	0,088+/-0,1 in / 0,22+/-0,02 cm
Mass	3,04 g

The tolerances are $\pm 5\%$ if another value is not indicated

Passing Elevator 4

Full length	8,36+/-0,25 in / 21,74+/-0,63 cm
Tip outer diameter	0,25 in / 0,635 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Tip length	2,7+/-0,06 in / 6,86+/-0,15 cm
	Tubing outer diameter	0,125 in / 0,32 cm
	Smaller tubing outer diameter	0,06 in / 0,15 cm
	Elevating tip height	1,62+/-0,03 in / 4,11+/-0,08 cm
	Elevating tip width	1,50+/-0,09 in / 3,81+/-0,23 cm
	Mass	3,56 g
The tolerances are $\pm 5\%$ if another value is not indicated		
Revision Accessory Kits, models: 97790-60, 97790-75, 97790-90 The Models 97790-60, 97790-75, and 97790-90 Revision Accessory Kits are intended for use when revising or repositioning an implanted Medtronic Neuromodulation Vectris SureScan MRI 1x8 percutaneous lead. With the exception of the lead, each kit contains all the sterile accessories supplied in the Vectris SureScan MRI lead kits according to lead length.		
Table 1. Package/Kit Configurations		

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Model Number	60 cm Stylets ^a	75 cm Stylets ^a	90 cm Stylets ^a	Guidewire	Passing Straw	Torque Wrench	Tunneling Rod/Tip	14-Gauge Needle – Straight, and Curved	Anchor Dispenser Tool with Bumpy Anchor	Anchor Handle with Bumpy Anchor	Anchor Removal Tool	Product literature
9779060	X			X	X	X	X	X	X	X	X	X
9779075		X		X	X	X	X	X	X	X	X	X
9779090			X	X	X	X	X	X	X	X	X	X

^a This includes 0.010 in (0.254 mm) straight and bent and 0.012 in (0.305 mm) straight and bent stylets.

Anchor with anchor dispenser tool

Full height	128,83 mm
Tip height	1,32 in / 3,35 cm
Tip outer diameter	0,188 in / 0,48 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Shaft outer diameter	0,35 in / 0,89 cm
	Handle outer diameter	0,45 in / 1,09 cm
	Mass of assy	10,13 g
	The tolerances are $\pm 5\%$ if another value is not indicated	
	Spare anchor on spare handle	
	Full height	2,42 in / 6,15 cm
	Tip height	1,32 in / 3,35 cm
	Tip outer diameter	0,188 in / 0,48 cm
	Shaft outer diameter	0,35 in / 0,89 cm
	Handle outer diameter	0,45 in / 1,09 cm
	Mass of assy	4,5 g
	The tolerances are $\pm 5\%$ if another value is not indicated	
	Anchor removal tool	
	Full length	1,62 $\pm$ 0,05 in / 4,11 $\pm$ 0,13 cm
	Full height	0,71 $\pm$ 0,05 in / 1,80 $\pm$ 0,13 cm
	Mass	1,33 g
	The tolerances are $\pm 5\%$ if another value is not indicated	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Guide wire

Length	28 cm
Outer diameter	0,84 mm
Mass	0,76 g

The tolerances are $\pm 5\%$ if another value is not indicated

Needle – curved tip, 14 gauge, 8.9 cm

Assy length	114,25 mm
Total length of needle	108,78 mm
Working length of needle	90,17 mm
Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,6 mm
Working length of mandrin	102,56 mm
Outer diameter of mandrin	1,39 mm
Hardness	≤ 92 HRC
Roughness	$Ra \leq 0,16 \mu m$
Bending strength of needle	$40 N \pm 0,1 N$
Blunting radius	$\leq 0,03 mm$
Mass	6,48 g

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Luer adaptor	Luer-lok
The tolerances are $\pm 5\%$ if another value is not indicated	
Needle – straight tip, 14 gauge, 8.9 cm	
Assy length	111,95 mm
Total length of needle	105,97 mm
Working length of needle	87,36 mm
Outer diameter of needle	1,950 mm – 2,150 mm
Minimum inner diameter of needle	1,6 mm
Working length of mandrin	99,44 mm
Outer diameter of mandrin	1,23 mm
Hardness	≤ 92 HRC
Roughness	$Ra \leq 0,16 \mu m$
Bending strength of needle	$40 N \pm 0,1 N$
Blunting radius	$\leq 0,03 mm$
Mass	6,38 g
Luer adaptor	Luer-lok
The tolerances are $\pm 5\%$ if another value is not indicated	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Straight Stylet, 0,25 mm

Dim L (wire length)	60 cm
	75 cm
	90 cm
Mass	1,6 g (60 cm)
	2 g (75 cm)
	2,4 g (90 cm)

The tolerances are $\pm 5\%$ if another value is not indicated

Bent Stylet 0,25 mm

Dim L (wire length)	60 cm
	75 cm
	90 cm
Mass	1,64 g (60 cm)
	2,05 g (75 cm)
	2,46 g (90 cm)

The tolerances are $\pm 5\%$ if another value is not indicated

Straight stylet 0,3 mm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Dim L (wire length)		60 cm
			75 cm
			90 cm
	Mass		1,78 g (60 cm)
			2,23 g (75 cm)
			2,68 g (90 cm)
	The tolerances are $\pm 5\%$ if another value is not indicated		
	Bent stylet 0,3 mm		
	Dim L (wire length)		60 cm
			75 cm
			90 cm
	Mass		1,82 g (60 cm)
			2,29 g (75 cm)
			2,75 g (90 cm)
	The tolerances are $\pm 5\%$ if another value is not indicated		
	Tunneling Rod		
	Tubing diameter		3,2 mm
	Outer diameter of mounting thread		3,11 mm
	Tubing length		8,14 in / 20,68 cm

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Mass	21,26 g
	Hardness	≤92 HRC
	Roughness	Ra≤3,2 μm
	The tolerances are ±5% if another value is not indicated	
	Tunneling Tube	
	Outer diameter	0,192+/-0,005 in / 4,9+/-0,1 mm
	Inner diameter	2,62 mm
	The tolerances are ±5% if another value is not indicated	
	Tunneling Tip	
	Diameter	3,92 mm
	Inner diameter of mounting thread	2,6 mm
	Tip length	1,01+/-0,06 in / 2,56 +/-0,15 cm
	Tip angle	15+/-3 °
	Blunting radius	≤0,03 mm
	Mass (with tunneling tube)	1,7 g
	Hardness	≤92 HRC
	Roughness	Ra≤3,2 μm
	The tolerances are ±5% if another value is not indicated	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Wrench, torque <table border="1"> <tr> <td>Hex tip outer diameter</td><td>From 0,40 in to 0,45 in / 0,10 to 0,11 cm</td></tr> <tr> <td>Screw outer diameter</td><td>0,434 in / 1,1 cm</td></tr> <tr> <td>Screw length</td><td>1,15+/-0,15 in / 2,92+/-0,38 cm</td></tr> <tr> <td>Hex tip length</td><td>,055+/-0,15 in / 1,40+/-0,38 cm</td></tr> <tr> <td>Mass</td><td>2,93 g</td></tr> <tr> <td>Hardness</td><td>48 HRC</td></tr> <tr> <td>Roughness</td><td>Ra≤3,2 μm</td></tr> </table> The tolerances are ±5% if another value is not indicated	Hex tip outer diameter	From 0,40 in to 0,45 in / 0,10 to 0,11 cm	Screw outer diameter	0,434 in / 1,1 cm	Screw length	1,15+/-0,15 in / 2,92+/-0,38 cm	Hex tip length	,055+/-0,15 in / 1,40+/-0,38 cm	Mass	2,93 g	Hardness	48 HRC	Roughness	Ra≤3,2 μm
Hex tip outer diameter	From 0,40 in to 0,45 in / 0,10 to 0,11 cm														
Screw outer diameter	0,434 in / 1,1 cm														
Screw length	1,15+/-0,15 in / 2,92+/-0,38 cm														
Hex tip length	,055+/-0,15 in / 1,40+/-0,38 cm														
Mass	2,93 g														
Hardness	48 HRC														
Roughness	Ra≤3,2 μm														
6. Description of the components (units) of medical device (if applicable)	I The Medtronic Models 3487A PISCES Quad Leads Package contents: - Lead Model 3487A – 1p. - Percutaneous extension – 1p. - Anchors: 2 wing – 1p., 3 wing – 1p., EZ – 1p., Twist lock – 1p. - Connector boot – 1p. - Guide wire – 1p. - Needle, 15 gauge (modified Tuohy) – 1p. - Screener cable – 1p. - Stylets: 0.25 mm bent at ≈ 20° angle – 1p., 0.25 mm straight – 1p., 0.25 mm short – 1p. - Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p. - Wrench, torque – 1p. - Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., indications														

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

and additional precautions – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 37082 Package Contents:

- Extension – 1p.
- Connector boots, 1 x 4 low-profile – 2p.
- Octapolar in-line neurostimulator plug – 1p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p., carrier – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

2. The Medtronic Model 37083 Package Contents:

- Extension – 1p.
- Connector boots, 1 x 4 low-profile – 2p.
- Octapolar in-line neurostimulator plug – 1p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p., carrier – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

3. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Flat Boot – 1p.
- 1x8 low profile closed – 1p.
- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

4. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

5. Model 3550-38 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Single-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product Literature: information letter (instructions for use) – 1p.

6. The Medtronic Model 3550-60 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Dual-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product literature: information letter (instructions for use) – 1p.

7. The Medtronic Model 3555-60 Package Contents:

- Passer 60 cm – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Tunneling obturator – 1p.
- Carrier – 1p.
- Product literature: information letter (instructions for use) – 1p.

8. Model 3550-43 Package Contents:

- 14-gauge 6" curved-tip needle – 1p.
- Product literature: information letter (instructions for use) – 1p.

II The Medtronic Models 3887 Quad Compact Leads Package contents:

- Lead Model 3887 – 1p.
- Percutaneous extension – 1p.
- Anchors: 2 wing – 1p., 3 wing – 1p., EZ – 1p., Twist lock – 1p.
- Connector boot – 1p.
- Guide wire – 1p.
- Needle, 15 gauge (modified Tuohy) – 1p.
- Screener cable – 1p.
- Stylets: 0.25 mm bent at $\approx 20^\circ$ angle – 1p., 0.25 mm straight – 1p., 0.25 mm short – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., indications and additional precautions – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

Accessories:

1. The Medtronic Model 37082 Package Contents:

- Extension – 1p.
- Connector boots, 1 x 4 low-profile – 2p.
- Octapolar in-line neurostimulator plug – 1p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p., carrier – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

2. The Medtronic Model 37083 Package Contents:

- Extension – 1p.
- Connector boots, 1 x 4 low-profile – 2p.
- Octapolar in-line neurostimulator plug – 1p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p., carrier – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

3. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.
- Flat Boot – 1p.
- 1x8 low profile closed – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

4. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

5. Model 3550-38 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Single-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product Literature: information letter (instructions for use) – 1p.

6. The Medtronic Model 3550-60 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Dual-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product literature: information letter (instructions for use) – 1p.

7. The Medtronic Model 3555-60 Package Contents:

- Passer 60 cm – 1p.
- Tunneling obturator – 1p.
- Carrier – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Product literature: information letter (instructions for use) – 1p.

8. Model 3550-43 Package Contents:

- 14-gauge 6" curved-tip needle – 1p.
- Product literature: information letter (instructions for use) – 1p.

III The Medtronic Models 3888 PISCES Quad Plus Leads Package contents:

- Lead Model 3888 – 1p.
- Percutaneous extension – 1p.
- Anchors: 2 wing – 1p., 3 wing – 1p., EZ – 1p., Twist lock – 1p.
- Connector boot – 1p.
- Guide wire – 1p.
- Needle, 15 gauge (modified Tuohy) – 1p.
- Screener cable – 1p.
- Stylets: 0.25 mm bent at $\approx 20^\circ$ angle – 1p., 0.25 mm straight – 1p., 0.25 mm short – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., indications and additional precautions – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 37082 Package Contents:

- Extension – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Connector boots, 1 x 4 low-profile – 2p.
- Octapolar in-line neurostimulator plug – 1p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p., carrier – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

2. The Medtronic Model 37083 Package Contents:

- Extension – 1p.
- Connector boots, 1 x 4 low-profile – 2p.
- Octapolar in-line neurostimulator plug – 1p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p., carrier – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

3. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.
- Flat Boot – 1p.
- 1x8 low profile closed – 1p.
- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

4. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

5. Model 3550-38 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Single-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product Literature: information letter (instructions for use) – 1p.

6. The Medtronic Model 3550-60 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Dual-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product literature: information letter (instructions for use) – 1p.

7. The Medtronic Model 3555-60 Package Contents:

- Passer 60 cm – 1p.
- Tunneling obturator – 1p.
- Carrier – 1p.
- Product literature: information letter (instructions for use) – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

8. Model 3550-43 Package Contents:

- 14-gauge 6" curved-tip needle – 1p.
- Product literature: information letter (instructions for use) – 1p.

IV The Medtronic Vectris SureScan MRI 1x8 Subcompact Models 977A160, 977A175, 977A190 Leads Package contents :

- Lead – 1p.
- Anchors: Injex – 1p., Injex with spare handle – 1p.
- Anchor removal tool – 1p.
- Guide wire – 1p.
- Needles (epidural, modified Tuohy): curved tip, 14 gauge, 8.9 cm – 1p., straight tip, 14 gauge, 8.9 cm – 1p.
- Stylets: bent, 0.30 mm diameter (inserted in the lead) – 1p., straight, 0.30 mm diameter – 1p., bent, 0.25 mm diameter – 1p., straight, 0.25 mm diameter – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Torque wrench – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.
- Flat Boot – 1p.
- 1x8 low profile closed – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

2. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

3. Model 3550-38 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Single-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product Literature: information letter (instructions for use) – 1p.

4. The Medtronic Model 3550-60 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Dual-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product literature: information letter (instructions for use) – 1p.

5. The Medtronic Model 3555-60 Package Contents:

- Passer 60 cm – 1p.
- Tunneling obturator – 1p.
- Carrier – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Product literature: information letter (instructions for use) – 1p.

6. The Medtronic Model 97791, 97792 Package Contents:

- Anchor with anchor dispenser tool – 2p.
- Anchor removal tool – 1p.
- Product literature: instructions for use – 1p.

7. The Medtronic Model 97790-60, 97790-75, 97790-90 Package Contents:

- Anchors: Injex – 1p., Injex with spare handle – 1p.
- Anchor removal tool – 1p.
- Guide wire – 1p.
- Needles (epidural, modified Tuohy): curved tip, 14 gauge, 8.9 cm – 1p., straight tip, 14 gauge, 8.9 cm – 1p.
- Stylets: bent, 0.30 mm diameter (inserted in the lead) – 1p., straight, 0.30 mm diameter – 1p., bent, 0.25 mm diameter – 1p., straight, 0.25 mm diameter – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Torque wrench – 1p.
- Product literature: information letter (instructions for use) – 1p.

8. Model 3550-43 Package Contents:

- 14-gauge 6" curved-tip needle – 1p.
- Product literature: information letter (instructions for use) – 1p.

9. The Medtronic Model 37081 Package Contents.:

- Extension – 1p.
- Connector boots, 1 x 8 low-profile – 2p.
- Setscrews – 2p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Tunneling tools: extension passer – 1p., obturator – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

V The Medtronic Vectris SureScan MRI 1x8 Compact Models 977A260, 977A275, 977A290 Leads Package contents:

- Lead – 1p.
- Injex Anchor – 1p.
- Anchor removal tool – 1p.
- Guide wire – 1p.
- Needles (epidural, modified Tuohy): curved tip, 14 gauge, 8.9 cm – 1p., straight tip, 14 gauge, 8.9 cm – 1p.
- Stylets: bent, 0.30 mm diameter (inserted in the lead) – 1p., straight, 0.30 mm diameter – 1p., bent, 0.25 mm diameter – 1p., straight, 0.25 mm diameter – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Torque wrench – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.
- Flat Boot – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- 1x8 low profile closed – 1p.
- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

2. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

3. Model 3550-38 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Single-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product Literature: information letter (instructions for use) – 1p.

4. The Medtronic Model 3550-60 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Dual-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product literature: information letter (instructions for use) – 1p.

5. The Medtronic Model 3555-60 Package Contents:

- Passer 60 cm – 1p.
- Tunneling obturator – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Carrier – 1p.
- Product literature: information letter (instructions for use) – 1p.

- 6. The Medtronic Model 97791, 97792 Package Contents:**
 - Anchor with anchor dispenser tool – 2p.
 - Anchor removal tool – 1p.
 - Product literature: instructions for use – 1p.

- 7. The Medtronic Model 97790-60, 97790-75, 97790-90 Package Contents:**
 - Anchors: Injex – 1p., Injex with spare handle – 1p.
 - Anchor removal tool – 1p.
 - Guide wire – 1p.
 - Needles (epidural, modified Tuohy): curved tip, 14 gauge, 8.9 cm – 1p., straight tip, 14 gauge, 8.9 cm – 1p.
 - Stylets: bent, 0.30 mm diameter (inserted in the lead) – 1p., straight, 0.30 mm diameter – 1p., bent, 0.25 mm diameter – 1p., straight, 0.25 mm diameter – 1p.
 - Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
 - Torque wrench – 1p.
 - Product literature: information letter (instructions for use) – 1p.

- 8. Model 3550-43 Package Contents:**
 - 14-gauge 6" curved-tip needle – 1p.
 - Product literature: information letter (instructions for use) – 1p.

- 9. The Medtronic Model 37081 Package Contents.:**
 - Extension – 1p.
 - Connector boots, 1 x 8 low-profile – 2p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

VI The Medtronic Vectris 1x8 Subcompact Models 977D160 and the Vectris 1x8 Compact Models 977D260 Leads **Package contents:**

- Lead – 1p.
- Guide wire – 1p.
- Needles (epidural, modified Tuohy): curved tip, 14 gauge, 8.9 cm – 1p., straight tip, 14 gauge, 8.9 cm – 1p.
- Stylets: bent, 0.30 mm diameter (inserted in the lead) – 1p., straight, 0.30 mm diameter – 1p., bent, 0.25 mm diameter – 1p., straight, 0.25 mm diameter – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

2. Model 3550-38 Package Contents:

- Passer – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Tunneling obturator – 1p.
- Single-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product Literature: information letter (instructions for use) – 1p.

4. The Medtronic Model 3550-60 Package Contents:

- Passer – 1p.
- Tunneling obturator – 1p.
- Dual-carrier obturator – 1p.
- 1x8 low-profile boot – 1p.
- Product literature: information letter (instructions for use) – 1p.

5. Model 3550-43 Package Contents:

- 14-gauge 6" curved-tip needle – 1p.
- Product literature: information letter (instructions for use) – 1p.

VII Medtronic Specify SureScan MRI Surgical Leads, Models 977C163, 977C190 Package contents:

- Lead – 1p.
- Anchors: two-wing – 2p., three-wing – 2p., – EZ – 2p.
- Passing elevator – 1p.
- Lead blank – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Torque wrench – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., list of

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 37081 Package Contents.:

- Extension – 1p.
- Connector boots, 1 x 8 low-profile – 2p.
- Setscrews – 2p.
- Tunneling tools: extension passer – 1p., obturator – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

2. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.
- Flat Boot – 1p.
- 1x8 low profile closed – 1p.
- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

3. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Product literature: instructions for use – 1p.

4. Model 3550-P4 Package Contents:

- Passing elevator 1 – 1p.
- Passing elevator 2 – 1p.
- Passing elevator 3 – 1p.
- Passing elevator 4 – 1p.
- Product Literature: instructions for use – 1p.

VIII Medtronic Specify SureScan MRI Surgical Leads, Models 977C265, and 977C290 Package contents:

- Lead – 1p.
- Anchors: two-wing – 2p., three-wing – 2p., – EZ – 2p.
- Passing elevator – 1p.
- Lead blank – 1p.
- Tunneling tools: tunneling rod – 1p., tunneling tip – 1p., tube – 1p.
- Torque wrench – 1p.
- Product literature: implantation manual – 1p., information for physician – 1p., indications information – 1p., list of patient stickers – 1p., MRI guidance – 1p.

Accessories:

1. The Medtronic Model 37081 Package Contents.:

- Extension – 1p.
- Connector boots, 1 x 8 low-profile – 2p.
- Setscrews – 2p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

- Tunneling tools: extension passer – 1p., obturator – 1p.
- Wrench, torque – 1p.
- Product literature: implantation manual – 1p., indications information – 1p.

2. The Medtronic Model 3550-24 Package Contents:

- 1x4 Boot – 2p.
- 1 x 8 Boot – 1p.
- 1 x 4 low profile Boot – 2p.
- 1 x 8 low profile Boot – 2p.
- Flat Boot – 1p.
- 1x8 low profile closed – 1p.
- Anchors: EZ – 1p., 2 wing – 1p., 3 wing – 1p.
- Product literature: guidance (instructions for use) – 1p.

3. The Medtronic Model 355531 Package Contents:

- Multi-lead trialing cable – 1p.
- Product literature: instructions for use – 1p.

4. Model 3550-P4 Package Contents:

- Passing elevator 1 – 1p.
- Passing elevator 2 – 1p.
- Passing elevator 3 – 1p.
- Passing elevator 4 – 1p.
- Product Literature: instructions for use – 1p.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

7. Description of all accessories, medical devices or non-medical devices, intended for use with the medical device submitted for registration (if applicable)	Neurostimulation Itrel 4 (Russian registration №P3H 2016/3601 from 25.01.2016) Itrel 4 shares the following system components and accessories currently approved with Itrel 3 and Restore family devices. SCS extensions and accessories - Model 37083 compatible with Itrel 4 Model 37704 and Restore family devices - Models 7489 compatible with Itrel 4 Model 37703, Itrel 3 SCS leads - Models 3986A, 3487A, 3887, 3888 External neurostimulator (ENS) - Model 37022 SureScan MRI system is designed to be compatible with the commercially available Medtronic neurostimulation system components that are approved to be used with the Restore system family, with the exception of requiring software updates to the external programming devices (Patient Recharger, Patient Programmer, and N'Vision Clinician Programmer (Russian registration № ΦC3 2011/11151 from 06.12.2011)) INS Models - PrimeAdvanced SureScan MRI Model 97702 - RestoreUltra SureScan MRI Model 97712 - RestoreAdvanced SureScan MRI Model 97713 - RestoreSensor SureScan MRI Model 97714 Leads - Vectris SureScan MRI 1x8 Subcompact leads Models 977A160, 977A175, and 977A190 - Vectris SureScan MRI 1x8 Compact Models 977A260, 977A275, and 977A290 - Vectris 1x8 Subcompact Model 977D160 trialing lead

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
Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

			European conformity – This symbol means that the device fully complies with European Directive AIMD 90/385/EEC	
			For USA audiences only	
			Official representative in EC	
			Manufacturer	
			Refer to IFU	
			Do not use if package is damaged or missing	
			Neurostimulation for chronic pain	
			Lead kit	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

			MR Conditional
			Package contents
			Lead
			Product documentation
			Accessories
			Sterile
9. List of materials of animal and (or) human origin, indicating		N/A	

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

information on their biocompatibility and safety, on selection of sources (donors), sampling, processing, storage and handling of these materials (if applicable)	
10. Information from the Russian State Pharma Products' Register about Pharma product contained in medical device (if applicable)	N/A
11. Sterilization method, Sterilization Validation, Packaging Validation for Sterile devices	Sterilization - Medtronic has sterilized the package contents according to the process indicated on the package label before shipment. This device is for single use only and is not intended to be resterilized. Do not Implant a component if the following circumstances have occurred: • The storage package has been pierced or altered because component sterility cannot be guaranteed and infection may occur. • The component shows signs of damage because the component may not function properly. • The use-by date has expired because component sterility cannot be guaranteed and infection may occur; also, neurostimulator battery longevity may be reduced and may require early replacement.

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	Sterilized by EO. Sterilization has been validated according to EN/ISO requirements. Annual validation at the contract facility is based upon EN 556-1/ISO 11135-1, Medical Devices Validation and routine control of ethylene oxide sterilization. Analysis of ethylene oxide residuals and its byproducts is in compliance with ISO 10993-7, Biological Evaluation of Medical Devices – part 7: Ethylene oxide sterilization residuals.
12. Requirements for maintenance and repair of a medical device	N/A
13. Procedure and conditions for the disposal or destruction of a medical device	Follow local regulations for proper disposal of this product. Do not dispose of this product in the unsorted municipal waste stream. This product contains materials that can harm the environment. See http://recycling.medtronic.com for instructions on proper disposal of this product.
14. Information about manufacturer of medical device (name, location address, mailing address, phone number, address of place of production)	1. Medtronic Neuromodulation, 800 53rd Ave. N.e., Minneapolis MN 55421, USA 2. Medtronic Neuromodulation, 7000 Central Ave. N.e., Minneapolis MN 55432, USA 3. Medtronic Puerto Rico Operations Co., Juncos, Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos, PR 00777, USA 4. Medtronic Europe Sarl, Route du Molliau 31, Case Postale, 1131 Tolochenaz, Switzerland 5. Medtronic Puerto Rico Operations Co., Villalba, Rd. 149, Km. 56.3, Call Box 6001, Villalba, PR 00766, USA
15. Class risk depending on the potential risk of a medical product, in accordance with the nomenclature	3

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

16. Type of contact with the body	3487A- Contact with the internal environment of the body > 30 days (Implantable device) 3887- Contact with the internal environment of the body > 30 days (Implantable device) 3888- Contact with the internal environment of the body > 30 days (Implantable device) 977A160- Contact with the internal environment of the body > 30 days (Implantable device) 977A175- Contact with the internal environment of the body > 30 days (Implantable device) 977A190- Contact with the internal environment of the body > 30 days (Implantable device) 977A260- Contact with the internal environment of the body > 30 days (Implantable device) 977A275- Contact with the internal environment of the body > 30 days (Implantable device) 977A290- Contact with the internal environment of the body > 30 days (Implantable device) 977D160- Contact with the internal environment of the body < 30 days (Temporary Implantable device) 977D260- Contact with the internal environment of the body < 30 days (Temporary Implantable device) 977C165- Contact with the internal environment of the body > 30 days (Implantable device) 977C190- Contact with the internal environment of the body > 30 days (Implantable device) 977C265- Contact with the internal environment of the body > 30 days (Implantable device) 977C290- Contact with the internal environment of the body > 30 days (Implantable device) Accessories: 37081- Contact with the internal environment of the body > 30 days (Implantable device) 37082- Contact with the internal environment of the body > 30 days (Implantable device) 37083- Contact with the internal environment of the body > 30 days (Implantable device) 3550-24- Contact with the internal environment of the body > 30 days (Implantable device) 3550-38- Contact with the internal environment of the body< 24 hours (Non-implanted Device) 3550-43- Contact with the internal environment of the body< 24 hours (Non-implanted Device) 3550-P4- Contact with the internal environment of the body< 24 hours (Non-implanted Device) 3550-60, 3555-60 - Contact with the internal environment of the body< 24 hours (Non-implanted Device) 355531- Contact with skin < 30 days (Non-implanted Device) 97790-60- Contact with the internal environment of the body< 24 hours (Non-implanted Device)
--	--

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	97790-75- Contact with the internal environment of the body< 24 hours (Non-implanted Device) 97790-90- Contact with the internal environment of the body< 24 hours (Non-implanted Device)
17. Methods and apparatus for disinfection and sterilization cleaning	N/A
18. Shelf life	3487A - 48 Month Shelf Life 3887 - 48 Month Shelf Life 3888 - 48 Month Shelf Life 977A160 - 48 Month Shelf Life 977A175 - 48 Month Shelf Life 977A190 - 48 Month Shelf Life 977A260 - 48 Month Shelf Life 977A275 - 48 Month Shelf Life 977A290 - 48 Month Shelf Life 977D160 - 48 Month Shelf Life 977D260 - 48 Month Shelf Life 977C165 - 48 Month Shelf Life 977C190 - 48 Month Shelf Life 977C265 - 48 Month Shelf Life 977C290 - 48 Month Shelf Life Accessories: 37081 - 48 Month Shelf Life 37082 - 48 Month Shelf Life 37083 - 48 Month Shelf Life 3550-24 - 48 Month Shelf Life 3550-38 -48 Month Shelf Life 3550-43 -48 Month Shelf Life

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	3550-P4 -48 Month Shelf Life 3555-60, 3550-60 -48 Month Shelf Life 355531 - 48 Month Shelf Life 97790-60 - 48 Month Shelf Life 97790-75 - 48 Month Shelf Life 97790-90 - 48 Month Shelf Life
19. Operating conditions	Single-Use, Do Not Reuse The devices can be used in normal conditions. Temperature range of operating room: 15°C~30°C. Temperature range for exploitation: from 32°C till 42°C, humidity from 80% till 100% (+-2%).
20. Conditions of transportation and storage	The device is transported by all modes of transport in covered vehicles in conformity with the rules of transportation of goods, operating in each mode. Store the leads between -34 °C (-30 °F and 57 °C (135 °F), humidity till 80%. Temperature transport -30°F to 135°F -34°C till 57°C, humidity till 80%. Temperatures outside of range above can damage components.
21. Installation requirements	N/A
22. Information provided by the manufacturer of the medical device, keys, access passwords, programs required for installation, adjustment, operation and maintenance of the	N/A

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

medical device	
23. Warranty	Limited Warranty- Areas outside the United States should contact their local Medtronic representative for exact terms of the Limited Warranty, In Russia please contact with Authorized Representative in Russia: Medtronic, LLC, 123317, Moscow, Presnenskaya Naberezhnaya., 10, Tel.: (495) 580-73-77,
24. The list of national and international regulatory documents/standards, which matches medical device	FDA (under Federal Food, Drug, and Cosmetic Act), EN ISO 13485 (Medical Devices – Quality Management Systems – Requirements for regulatory purposes), EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals EN ISO 11135-1 Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices EN ISO 11607-1 Packaging for terminally sterilized medical Devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems EN ISO 11607-2 Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes EN ISO 11737-2 Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process EN ISO 14971 Medical devices - Application of risk management to medical devices EN1041 Information Supplied By The Manufacturer With Medical Devices EN45502-1 Active Implantable Medical Devices – Part 1: General Requirements For Safety, Marking and information To Be Provided By The Manufacturer EN556-1 Sterilization of medical devices – Requirements for medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilized medical devices EN62366 Medical devices — Application of usability engineering to medical devices

Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com

	EN980 Graphical Symbols For Use In The Labeling Of Medical Devices
25. Customer Service	For all questions and solving of problems with use of the products Medtronic in Russia please contact with Authorized Representative in Russia: Medtronic, LLC, 123317, Moscow, Presnenskaya Naberezhnaya., 10 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru

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Medtronic

Medtronic Pain Therapy
Using neurostimulation for chronic pain


Leczenie bólu przy użyciu urządzeń firmy Medtronic
Stosowanie neurostymulacji przy bólu przewlekłym
Léčba bolesti od společnosti Medtronic
Použití neurostimulace pro léčbu chronické bolesti
Medtronic fájdalomkezelés
Neurostimuláció alkalmazása a krónikus fájdalom kezelésében
Обезболивающая терапия Medtronic
Использование нейростимуляции при хронической боли
Liečba bolesti od spoločnosti Medtronic
Používanie neurostimulácie na liečbu chronickej bolesti
Medtronic Ağrı Tedavisi
Kronik ağrı için nörostimülasyon kullanımı
Soluții Medtronic pentru terapia durerii
Utilizarea neurostimulării pentru durerile cronice

Information for prescribers • Informacje dla lekarzy zlecających • Informace
pro lékaře předepisující léčbu • Tájékoztató az előíró orvosnak •
Информация для врачей • Informácie pre predpisujúcich lekárov • Reçete
yazanlar için bilgiler • Informații pentru medicii prescriptori

USA Rx only

CE0123

*As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic, Inc.*

Ashley M. Luloff
November 15, 2019.

Medtronic and Medtronic logo are trademarks of Medtronic. Irit[™] is a trademark of a Medtronic company.

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Table of contents

Contraindications	5
Warnings	5
Precautions	8
Use in specific populations	8
Physician training	8
Storage and sterilization	8
System implant	8
Clinician programming	11
Patient programming and patient control devices	12
Charging system (for rechargeable neurostimulators)	12
Patient activities	13
Component disposal	14
Individualization of Treatment	14
Adverse events summary	15
Patient Counseling Information	16

Refer to the Indications sheet for indications and related information.

Refer to the MRI Guidelines for Medtronic Neurostimulation Systems for Chronic Pain Instructions for use manual for the MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan.

Refer to the device implant manual for device description, package contents, device specifications, and instructions for use.

Refer to the System Eligibility, Battery Longevity, Specifications reference manual for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.

USE Refer to the clinical summary booklet for information on the clinical study results of the neurostimulation system and individualization of treatment.

Appendix A: Electromagnetic Interference 17

Contraindications 17

Warnings 17

Precautions 19

Notes 21

4 English Information for prescribers 2018-07-31

Contraindications

Diathermy - Do not use shortwave diathermy, microwave diathermy or therapeutic ultrasound diathermy (all now referred to as diathermy) on patients implanted with a neurostimulation system. Energy from diathermy can be transferred through the implanted system and can cause tissue damage at the location of the implanted electrodes, resulting in severe injury or death. Refer to Appendix A: Electromagnetic Interference for further information.

Warnings

Electromagnetic Interference (EMI) - Electromagnetic interference is a field of energy generated by equipment found in the home, work, medical or public environments that is strong enough to interfere with neurostimulator function. Neurostimulators include features that provide protection from electromagnetic interference. Most electrical devices and magnets encountered in a normal day are unlikely to affect the operation of a neurostimulator. However, sources of strong electromagnetic interference can result in the following:

- Serious patient injury or death, resulting from heating of the implanted components of the neurostimulation system and damage to surrounding tissue.
- System damage, resulting in a loss of or change in symptom control and requiring surgical replacement.
- Operational changes to the neurostimulator, causing it to turn on or off (particularly in neurostimulators enabled for magnet use), or to reset to power-on-reset (POR) settings, resulting in loss of stimulation, return of symptoms, and in the case of POR, requiring reprogramming by a clinician.
- Unexpected changes in stimulation, causing a momentary increase in stimulation or intermittent stimulation, which some patients have described as a jolting or shocking sensation. Although the unexpected change in stimulation may feel uncomfortable, it does not damage the device or injure the patient directly. In rare cases, as a result of the unexpected change in stimulation, patients have fallen down and been injured.

Refer to Table 1, starting on page 8, and Appendix A: Electromagnetic Interference for information on sources of EMI, the effect of EMI on the patient and the neurostimulation system, and instructions on how to reduce the risk from EMI.

For information about the effects of EMI on programming, refer to "Telemetry signal disruption from EMI" on page 11.

Table 1. Potential effects of EMI from devices or procedures

Device or procedure	Serious patient injury	Device damage	Device turns off/on	Momentary increase in stimulation	Intermittent stimulation	For guidelines
Bone growth stimulators		X		X	X	page 19
CT scans				X		page 17
Defibrillation/ cardioversion	X	X		X	X	page 17
Dental drills and ultrasonic probes		X				page 19
Diathermy, therapeutic	X	X			X	page 17
Electrocautery	X	X				page 18
Electrolysis		X			X	page 19
Electromagnetic field devices: (eg, arc welding, power stations)			X	X	X	page 20
High-output ultrasonics / lithotripsy		X				page 18
Household items			X	X		page 21
Laser procedures		X				page 20
Magnetic resonance imaging (MRI)	X	X	X	X	X	page 7, page 18
Psychotherapeutic procedures		X	X	X	X	page 20
Radiation therapy		X				page 20
Radio-frequency (RF) / microwave ablation	X	X			X	page 18
Theft detector			X	X	X	page 18

Table 1. Potential effects of EMI from devices or procedures (continued)

Device or procedure	Serious patient injury	Device damage	Device turns off/on	Momentary increase in stimulation	Intermittent stimulation	For guidelines
Therapeutic magnets			X			page 21
Therapeutic ultrasound	X	X			X	page 17
Transcutaneous electrical nerve stimulation (TENS)			X	X		page 20

Magnetic resonance imaging (MRI) - Refer to the MRI guidelines manual associated with this product for the MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan. Do not conduct an MRI examination on a patient with any implanted neurostimulation system component until you read and fully understand all MRI information in the MRI guidelines manual. MR scans performed under different conditions can result in patient injury or damage to the implantable device.

Case damage - If the neurostimulator case is ruptured or pierced due to outside forces, severe burns could result from exposure to the battery chemicals.

Neurostimulator interaction with implanted cardiac devices - When a patient's medical condition requires both a neurostimulator and an implanted cardiac device (eg, pacemaker, defibrillator), physicians involved with both devices (eg, neurologist, neurosurgeon, cardiologist, cardiac surgeon) should discuss the possible interactions between the devices before surgery. To minimize or prevent the effects described below, implant the devices on opposite sides of the body and follow any additional instructions.

- Defibrillation therapy from an implanted defibrillator may damage the neurostimulator.
- The electrical pulses from the neurostimulation system may interact with the sensing operation from a cardiac device and could result in an inappropriate response of the cardiac device. To minimize or prevent the cardiac device from sensing the neurostimulator output, program the neurostimulator to a bipolar configuration and to a minimum rate of 60 Hertz. Program the cardiac device to bipolar sensing.

See also: "Programmer interaction with other active implanted devices" on page 11.

Precautions

Use in specific populations

Use in specific populations - The safety and effectiveness of this therapy has not been established for:

- pregnancy, unborn fetus, or delivery.
- pediatric use (patients under the age of 18).

Physician training

Implanting physicians - Implanting physicians should be experienced in spinal procedures and should review the procedures described in the implant manual before surgery.

Prescribing physicians - Prescribing physicians should be experienced in the diagnosis and treatment of chronic intractable pain and should be familiar with using the neurostimulation system.

Storage and sterilization

Component packaging - Do not implant a component if the following circumstances have occurred:

- The storage package has been pierced or altered because component sterility cannot be guaranteed and infection may occur.
- The component shows signs of damage because the component may not function properly.
- The use-by date has expired because component sterility cannot be guaranteed and infection may occur; also, neurostimulator battery longevity may be reduced and may require early replacement.

Sterilization - Medtronic has sterilized the package contents according to the process indicated on the package label before shipment. This device is for single use only and is not intended to be re-sterilized.

Storage temperature: leads and extensions - Do not store or transport the lead or extension above 57 °C (135 °F) or below -34 °C (-30° F). Temperatures outside this range can damage components.

Storage temperature: neurostimulators - Do not store or transport the neurostimulator above 52 °C (125 °F) or below -18 °C (0 °F). Temperatures outside this range can damage components.

System implant

Compatibility, all components - Follow these guidelines when selecting system components:

- **Medtronic components:** For proper therapy, use only Medtronic neurostimulation components that are compatible or specified in an intended use statement (if

present).

Components are compatible when the following conditions are met:

- Components have the same indication.
- For implanted components, the contact spacing and the number of electrode contacts at the connections for the lead and extension/neurostimulator or extension and neurostimulator are the same.

For each product, refer to the indication insert(s) and the shipping label artwork for this information.

- **Non-Medtronic components:** No claims of safety or efficacy are made with regard to the compatibility of using components that are not supplied by Medtronic. Refer to the non-Medtronic documentation for information.

Compatibility, bifurcated extension - When using a bifurcated extension with a different number of contacts than the neurostimulator, select leads of similar impedance. Using a low impedance lead with a standard impedance lead may result in perceived stimulation differences. The effect of dissimilar leads is unknown. (A low impedance lead [impedance below 2 ohms / 10 cm] is a Medtronic product name or model number designated by "Z," or "LZ," or "low impedance.")

Component handling - Handle the implantable components of this system with extreme care. These components may be damaged by excessive traction or sharp instruments, which may result in intermittent or loss of stimulation, requiring surgical replacement. Refer to the appropriate implant manual for additional instructions.

Extension routing for multiple leads - When multiple leads are implanted, route the lead-extensions so the area between them is minimized (Figure 1). If the lead-extensions are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (tolling or shocking sensation).


Figure 1. Routing for multiple lead-extensions.

Neurostimulator location - Select a location that is:

- a minimum of 20 cm (8 in) away from another neurostimulator to minimize telemetry interference and possible inappropriate therapy.
- on the opposite side of the body from another active implanted device (eg, pacemaker, defibrillator) to minimize possible interaction between the devices.
- away from bony structures (eg, 3 – 4 cm [1.2 – 1.6 in]) to minimize discomfort at the neurostimulator site.
- away from areas of restriction or pressure to minimize the potential for skin erosion and patient discomfort.
- in an area accessible to the patient for proper operation of a patient control device (ie, patient programmer, control magnet, radio-frequency transmitter).

Device inversion - To minimize the risk of device inversion or rotation:

- Do not make the neurostimulator pocket any larger than necessary to fit the neurostimulator and excess lead or extension.
- Use the suture holes in the connector block to secure the neurostimulator to the muscle fascia per the implant manual.

Device inversion may also occur from erosion of subcutaneous tissue within the pocket, pocket fluid, or hematoma. See also "Patient activities" sections: "Activities requiring excessive bending or stretching" and "Component manipulation by patient."

Refer to the appropriate implant manual for additional instructions.

Patient detoxification - Before placing leads for test stimulation, patients should be detoxified from narcotics. If patients are not detoxified, test stimulation may not be properly assessed.

Clinician programming

Amplitude lower limit - During programming, ensure that the lower limit for amplitude is set to the lowest value available. If the patient cannot set the amplitude to the lowest value available before turning off the neurostimulator, the patient may experience unexpected or uncomfortable stimulation.

Parameter adjustment - To prevent possible uncomfortable or unexpected stimulation (ticking or shocking sensation):

- program parameter changes in small increments above the perception threshold (the amplitude value(s) at which the patient first perceives paresthesia).
- decrease the amplitude(s) to 0.5 V before:
 - changing groups.
 - changing electrode polarities.
 - connecting or disconnecting the cable to the screener or external neurostimulator.
 - replacing the external neurostimulator or screener batteries.
 - turning on the neurostimulator or screener.

Programmer interaction with a cochlear implant - When the patient has a cochlear implant, minimize or eliminate the potential for unintended audible clicks during telemetry by keeping the external portion of the cochlear system as far from the programming components as possible, or by turning off the cochlear implant during programming.

Programmer interaction with flammable atmospheres - The programmer is not certified for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. The consequences of using the programmer near flammable atmospheres are unknown.

Programmer interaction with other active implanted devices - When a patient has a neurostimulator and another active implanted device (eg, pacemaker, defibrillator, neurostimulator):

- the radio-frequency (RF) signal used to program these devices may reset or reprogram the other device.
- the magnet in a cardiac programmer may activate magnetically controlled functions of some neurostimulators.

To verify that inadvertent programming did not occur, clinicians familiar with each device should check the programmed parameters of each device before the patient is discharged from the hospital and after each programming session of either device (or as soon as possible after these times).

Also, inform patients to contact their physician immediately if they experience symptoms that could be related to either device or to the medical condition treated by either device.

Test stimulation with 2 screener devices or external neurostimulators - When programming test stimulation on a patient who has 2 leads implanted, do not operate 2 screener devices or external neurostimulators simultaneously. The signals from simultaneously operated screener devices and external neurostimulators can interfere with each other and result in incorrect results.

Telemetry signal disruption from EMI - Do not attempt telemetry near equipment that may generate electromagnetic interference (EMI). If EMI disrupts programming, move the

programmer away from the likely source of EMI. Examples of sources of EMI are magnetic resonance imaging (MRI), lithotripsy, computer monitors, cellular telephones, x-ray equipment, and other monitoring equipment.

Patient programming and patient control devices

Group selection - Patients should select the group recommended by the clinician for the current activity or posture. Use of another group may result in uncomfortable or unexpected stimulation (jolting or shocking) when stimulation is turned on.

Neurostimulator battery depletion - For the Irit 3 Neurostimulator, patients with a very low perception threshold (the amplitude at which the patient first perceives paresthesia) may feel the stimulation intensity fluctuate as the battery nears depletion. To compensate for this fluctuation, patients may need to decrease or increase the amplitude more often to maintain the desired level of symptom control.

Patient control devices may affect other implanted devices - Do not place the patient control device (ie, patient programmer, control magnet, radio-frequency transmitter) over another active implanted medical device (eg, pacemaker, defibrillator, another neurostimulator). The patient control device could unintentionally change the operation of the other device.

Patient magnet disable - If the neurostimulator magnet control has been programmed to off, the patient must have a patient programmer to turn the neurostimulator on and off.

Patient magnet may damage items - Do not place the patient magnet on or near computer monitors, magnetic storage disks or tapes, televisions, credit cards, or other items affected by strong magnetic fields. If the patient magnet is too close, these items may be damaged.

Patient device handling - To avoid damaging the device, do not immerse it in liquid; do not clean it with bleach, nail polish remover, mineral oil, or similar substances; and do not drop it or mishandle it in a way that may damage it.

Patient device use - When operating an external neurostimulator, patient programmer, or charging system, use special care near flammable or explosive atmospheres. An interaction between the flammable or explosive atmospheres and the battery in the device could occur. The consequences of using a battery-powered device near flammable or explosive atmospheres are unknown.

Screener or radio-frequency (RF) transmitter - Advise the patient not to adjust the physician controls. If the patient adjusts these controls, stimulation could be changed, resulting in uncomfortable or loss of stimulation.

Charging system (for rechargeable neurostimulators)

Wound contact - DO NOT use the recharger on an unhealed wound. The recharger, antenna, and belt are not sterile, and contact with the wound may cause an infection.

Recharger use - Check for skin irritation or redness near the neurostimulator during recharging. Do not sit or lie on the antenna, or apply excessive pressure to the antenna. Take periodic breaks during prolonged recharging. Although no direct cause and effect has been established, some patients have reported heating sensation, discomfort, blistering not caused by heating, skin irritation, or redness near the implanted

neurostimulator during or after recharging. Contributing factors may include excessive pressure on the antenna, prolonged recharging periods, or individual patient physiological factors.

Low battery charge level - Advise patients to charge the neurostimulator when a Low battery () screen is displayed on the patient programmer or recharger; this prevents the battery from overdischarging. If the neurostimulator battery is allowed to overdischarge, the patient cannot charge the neurostimulator; however, the clinician may be able to restore the battery function using the Physician Recharge Mode on the recharger (refer to the troubleshooting section of the software manual).

Allowing the neurostimulator battery to overdischarge will permanently affect the neurostimulator in one of the following ways:

- Battery function is restored; however, charging sessions may be more frequent because battery capacity has been reduced.
- Battery function is not restored and the neurostimulator must be surgically replaced. Battery function is not restored when:
 - the neurostimulator battery is permanently damaged.
 - the neurostimulator battery has been overdischarged and restored twice before. The third time the battery is overdischarged, the neurostimulator will reach end of service. Surgery is required to replace the neurostimulator.

Patient activities

Activities requiring excessive twisting or stretching - Patients should avoid activities that may put undue stress on the implanted components of the neurostimulation system.

Activities that include sudden, excessive, or repetitive bending, twisting, bouncing, or stretching can cause component fracture or dislodgement. Component fracture or dislodgement may result in loss of stimulation, intermittent stimulation, stimulation at the fracture site, and additional surgery to replace or reposition the component. Spinal cord stimulation patients, in particular, should avoid excessive bending of the torso.

Antenna attachment - Do not use the adhesive disc to attach the antenna over the neurostimulator incision until the incision heals. The adhesive discs may irritate or reopen the incision.

Antenna use - If swelling or redness occurs near the adhesive disc, advise the patient to contact the clinician before using the antenna again. Swelling or redness may indicate an infection or an allergic reaction to the adhesive disc or antenna.

Component manipulation by patient - Patients should avoid manipulating or rubbing the neurostimulation system through the skin. Manipulation may cause component damage, lead dislodgement, skin erosion, or stimulation at the implant site. Manipulation may also cause device inversion, making a rechargeable neurostimulator impossible to charge.

Scuba diving or hyperbaric chambers - Patients should not dive below 10 meters (33 feet) of water or enter hyperbaric chambers above 2.0 atmospheres absolute (ATA). Pressures below 10 meters (33 feet) of water (or above 2.0 ATA) could damage the neurostimulation system. Before diving or using a hyperbaric chamber, patients should discuss the effects of high pressure with their physician.

Skydiving, skiing, or hiking in the mountains - High altitudes should not affect the neurostimulator, however, the patient should consider the movements involved in any planned activity and take precaution to avoid putting undue stress on the implanted system.

Patients should be aware that during skydiving, the sudden jerking that occurs when the parachute opens may cause lead dislodgement or fractures, which may require surgery to repair or replace the lead.

Unexpected changes in stimulation - Electromagnetic interference, postural changes, and other activities may cause a perceived increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation); therefore, patients should reduce the amplitude to the lowest setting and turn off the neurostimulator before engaging in activities that could be unsafe for themselves or others if they received an unexpected jolt or shock (eg, driving, operating power tools). Patients should discuss these activities with their physician.

Component disposal

When explanting a device (eg, replacement, cessation of therapy, or postmortem), or when disposing of accessories, follow these guidelines:

- If possible, return the explanted device with completed paperwork to Medtronic for analysis and disposal. See the back cover for the mailing addresses.
- To allow for device analysis, do not autoclave the device or expose the device to ultrasonic cleaners.
- Dispose of any unreturned components according to local environmental regulations; in some countries, explanting a battery-powered implantable device is mandatory.
 - Do not incinerate or cremate the neurostimulator because it may explode if subjected to these temperatures.
 - Do not reuse any implantable device or implantable accessory after exposure to body tissues or fluids because the functionality of the component cannot be guaranteed.

Individualization of Treatment

Best results are achieved when the patient is fully informed about the therapy risks and benefits, surgical procedure, follow-up requirements, and self-care responsibilities. Maximum benefits from the neurostimulation system require long-term postsurgical management.

Patient selection

Select patients carefully to ensure that:

- their symptoms are of physiological origin.
- they are appropriate candidates for surgery.
- they can properly operate the system.
- they received satisfactory results from a trial or test stimulation period.

Use in specific populations

Use in specific populations - The safety and effectiveness of this therapy has not been established for:

- pregnancy, unborn fetus, or delivery.
- pediatric use (patients under the age of 18).

Long-term effectiveness of neurostimulation

The long-term effectiveness of neurostimulation has been documented. Long-term clinical data regarding the efficacy of Medtronic neurostimulation systems is not yet available. Not all patients realize long-term benefits from neurostimulation.

Adverse events summary

The implantation of a spinal cord stimulation system involves risks that are similar to other spinal procedures. In addition to those risks associated with surgery, the following adverse events may occur with implantation or use of a neurostimulation system:

- Allergic or immune system response to the implanted materials
- Infection
- Lead, extension, or neurostimulator erosion through the skin or migration
- Leakage of cerebrospinal fluid
- Loss of pain relief may return patients to their underlying pain condition
- Patients on anticoagulation therapies may be at greater risk for postoperative complications such as hematomas that can result in paralysis
- Persistent pain at the neurostimulator site
- Placement of the epidural lead-extension is a surgical procedure that may expose patients to risks of epidural hemorrhage, hematoma, or paralysis
- Radicular chest wall stimulation
- Seroma or hematomas at the neurostimulator site
- Change in stimulation, possibly related to cellular changes around the electrode(s), shifts in electrode position, loose electrical connections, lead or extension fractures, which has been described by some patients as uncomfortable stimulation (jolting or shocking sensation).
- Formation of reactive tissue around the lead in the epidural space can result in delayed spinal cord compression and paralysis, requiring surgical intervention. Time to onset can range from weeks to many years after implant.
- Stimulation-dependent gastrointestinal symptoms such as nausea, diarrhea, incontinence, or constipation
- Stimulation-dependent bladder symptoms such as urinary retention, incontinence, or frequency

USA For a comprehensive summary of adverse events, refer to the clinical summary.

Patient Counseling Information

Physicians should provide patients and caregivers with information about the following:

- The components of the neurostimulation system: lead, extension, and neurostimulator
- Instructions for using the neurostimulation system
- The indications, contraindications, warnings, precautions, and adverse events for a neurostimulation system

Physicians should also instruct patients to do the following:

- Always inform any health care personnel that they have an implanted neurostimulation system before any procedure is begun.
- Contact their physician if they notice any unusual symptoms or signs.

Appendix A: Electromagnetic interference

Please review Electromagnetic Interference (EMI) under "Warnings" on page 5 and Table 1: "Potential effects of EMI from devices or procedures" on page 8.

Before any medical procedure is begun, patients should always inform any health care personnel that they have an implanted neurostimulation system. The potential for the following effects results from an interaction of the neurostimulation system and equipment—even when both are working properly.

Contraindications

Diathermy - Do not use shortwave diathermy, microwave diathermy or therapeutic ultrasound diathermy (all now referred to as diathermy) on patients implanted with a neurostimulation system. Energy from diathermy can be transferred through the implanted system and can cause tissue damage at the location of the implanted electrodes, resulting in severe injury or death.

Diathermy can also damage the neurostimulation system components, resulting in loss of therapy and requiring additional surgery for system explantation and replacement. Advise your patient to inform all their health care professionals that they should not be exposed to diathermy treatment.

Injury to the patient or damage to the device can occur during diathermy treatment when:

- the neurostimulation system is turned on or off.
- diathermy is used anywhere on the body—not just at the location of the neurostimulation system.
- diathermy delivers heat or no heat.
- any component of the neurostimulation system (lead, extension, neurostimulator) remains in the body.

Warnings

EMI from the following medical procedures or equipment may damage the device, interfere with device operation, or cause harm to the patient. If these procedures are required, follow the guidelines below:

CT scans - Prior to the patient undergoing a CT scan, program the neurostimulator to 0 V and turn the neurostimulator off. If these guidelines are not followed, the patient may experience a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).

Defibrillation or cardioversion - When a patient is in ventricular or atrial fibrillation, the first consideration is patient survival. External defibrillation or cardioversion can damage a neurostimulation system and cause induced currents in the lead-extension portion of the neurostimulation system that can injure the patient. Minimize the current flowing through the neurostimulation system by following these guidelines:

- Position defibrillation paddles as far from the neurostimulator as possible.
- Position defibrillation paddles perpendicular to the neurostimulation system.
- Use the lowest clinically appropriate energy output (well seconds).

After defibrillation, confirm the neurostimulation system is functioning as intended.

Electrocautery - If electrocautery is used near an implantable device or contacts a device or insertion needle, the following effects may occur:

- The tissue surrounding the insertion-needle (during placement of a percutaneous lead) may be damaged.
- The insulation on the lead or extension may be damaged, resulting in component failure or induced currents into the patient that may damage tissue or stimulate or shock the patient.
- The neurostimulator may be damaged, output may be temporarily suppressed or increased, or stimulation may stop because parameters were changed to power-on-reset settings (eg. output off, amplitude 0.0 V).

When electrocautery is necessary, follow these precautions:

- Before using electrocautery, turn off the neurostimulator.
- Disconnect any cable connecting the lead or extension to a screener or external neurostimulator.
- Use only bipolar cautery.
- If unipolar cautery is necessary:
 - use only a low-voltage mode.
 - use the lowest possible power setting.
 - keep the current path (ground plate) as far from the neurostimulator, extension, and lead as possible.
 - do not use full-length operating room table grounding pads.
- After using electrocautery, confirm that the neurostimulator is functioning as intended.

High-output ultrasonics or lithotripsy - Use of high-output ultrasonic devices, such as electrohydraulic lithotriptors, is not recommended for patients who have an implanted neurostimulation system. If lithotripsy must be used, do not focus the beam within 15 cm (6 in) of the neurostimulator.

Magnetic resonance imaging (MRI) - Refer to the MRI guidelines manual associated with this product for the MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan. Do not conduct an MRI examination on a patient with any implanted neurostimulation system component until you read and fully understand all MRI information in the MRI guidelines manual. MR scans performed under different conditions can result in patient injury or damage to the implantable device.

Radio-frequency or microwave ablation - Safety has not been established for radio-frequency (RF) or microwave ablation in patients who have an implanted neurostimulation system. Induced electrical currents may cause heating, especially at the lead electrode site, resulting in tissue damage.

Theft detectors and security screening devices - Advise patients to use care when approaching theft detector and security screening devices (such as those found in airports, libraries, and some department stores). When approaching these devices, patients should do the following:

1. If possible, patients should request to bypass these devices. Patients should show the security personnel their patient identification card for the neurostimulator and

request a manual search. Security personnel may use a handheld security wand but patients should ask the security personnel not to hold the security wand near the neurostimulator any longer than is absolutely necessary. Patients may wish to ask for another form of personal search.

2. If patients must pass through the theft detector or security screening device, they should turn off their neurostimulator, approach the center of the device and walk through normally (Figure 2).
 - a. If two security gates are present, they should walk through the middle, keeping as far from each gate as possible.
 - b. If one gate is present, they should walk as far from it as possible.

Note: Some theft detectors may not be visible.

3. Patients should proceed through the security screening device. They should not linger near or lean on the security screening device.


Figure 2. Approaching security gates.

4. After patients pass through the security screening device, they should turn on their neurostimulator.

Precautions

EMI from the following equipment is unlikely to affect the neurostimulation system if the guidelines below are followed:

Bone growth stimulators - Keep external magnetic field bone growth stimulator coils 45 cm (18 in) away from the neurostimulation system. When using either an implantable or external bone growth stimulator, ensure that both the bone stimulator and neurostimulator are working as intended.

Dental drills and ultrasonic probes - Turn off the neurostimulator. Keep the drill or probe 15 cm (6 in) away from the neurostimulator.

Electrolysis - Turn off the neurostimulator. Keep the electrolysis wand at least 15 cm (6 in) away from the neurostimulator.

Electromagnetic field devices - Patients should exercise care or avoid the following equipment or environments:

- Antenna of citizens band (CB) radio or ham radio
- Electric arc welding equipment
- Electric induction heaters used in industry to bend plastic
- Electric steel furnaces
- High-power amateur transmitters
- High-voltage areas (safe if outside the fenced area)
- Linear power amplifiers
- Magnetic degaussing equipment
- Magnets or other equipment that generates strong magnetic fields
- Microwave communication transmitters (safe if outside the fenced area)
- Perfusion systems
- Resistance welders
- Television and radio transmitting towers (safe if outside the fenced area)

If patients suspect that equipment is interfering with neurostimulator function, they should do the following:

1. Move away from the equipment or object.
2. If possible, turn off the equipment or object.
3. Then, if necessary, use the control magnet or patient programmer to return the neurostimulator to the desired on and off state.
4. Inform the equipment owner or operator of the occurrence.

If the above actions do not resolve the effects of the interference, or the patients suspect that their therapy is not effective after exposure to EMI, they should contact their physician.

Laser procedures - Turn off the neurostimulator. Keep the laser directed away from the neurostimulation system.

Psychotherapeutic procedures - Safety has not been established for psychotherapeutic procedures using equipment that generates electromagnetic interference (eg, electroconvulsive therapy, transcranial magnetic stimulation) in patients who have an implanted neurostimulation system. Induced electrical currents may cause heating, especially at the lead electrode site, resulting in tissue damage.

Radiation therapy - High-radiation sources should not be directed at the neurostimulator. High-radiation exposure may temporarily interfere with neurostimulator operation and may damage the neurostimulator. Damage may not be immediately apparent. To limit device exposure, use appropriate shielding or other measures, such as making beam angle adjustments to avoid the device.

Transcutaneous electrical nerve stimulation - Do not place transcutaneous electrical nerve stimulation (TENS) electrodes so that the TENS current passes over any part of the neurostimulation system. If patients feel that the TENS may be interfering with the

implanted neurostimulator, patients should discontinue using the TENS until they talk with their doctor.

Notes

Household Items - Most household appliances and equipment that are working properly and grounded properly will not interfere with the neurostimulation system. The following equipment is generally safe if patients follow these guidelines:

- **Freezer, refrigerator, or storm door magnets that hold the door closed:** Do not lean against the magnetic strip of the door.
- **Radio-frequency sources (AM/FM radios, analog and digital cellular telephones, cordless and conventional telephones):** Keep these items at least 10 cm (4 in) away from the implanted neurostimulator.
- **Stereo speakers and radios for the home or car:** Do not sit or carry speakers or radios near the neurostimulator.
- **Sewing machines or salon hair dryer:** Keep the neurostimulator away from the motors.
- **Computer disk drives:** Keep the neurostimulator away from the disk drives.
- **Induction ranges:** Keep the neurostimulator away from the burners while the burners are turned on.
- **Power tools:** Keep the motor away from the neurostimulator, lead, and extension.

Other medical procedures - EMI from the following medical procedures is unlikely to affect the neurostimulation system:

- **Diagnostic ultrasound (eg, carotid scan, doppler studies)**

Note: To minimize potential image distortion, turn off the neurostimulator and keep the transducer 15 cm (6 in) away from the neurostimulation system.

- **Diagnostic x-rays or fluoroscopy**

Notes:

- To minimize potential image distortion, turn off the neurostimulator.
- Tight pressure such as used during mammography may damage the neurostimulator or disconnect the neurostimulation system components, which may require surgery to reconnect or replace components. During x-ray procedures that require external compression around implanted components, the x-ray equipment should be adjusted to limit the amount of pressure exerted on the neurostimulator.

- **Magnetoencephalography (MEG)**

- **Positron emission tomography (PET) scans**

Therapeutic magnets (eg, magnetic mattresses, blankets, wrist wraps, elbow wraps) - Keep the magnet at least 25 cm (10 in) away from the neurostimulator. Magnetic fields of 10 gauss or less will generally not affect the neurostimulator.

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
As a registered Notary Public for the State of Minnesota,
 I hereby verify that this is a true and accurate copy of the
 original document. The original is filed at Medtronic, Inc.

Ashley M. Luloff
 November 15, 2019.

Indications

Neurostimulation for spinal cord stimulation (SCS) - The Medtronic SCS neurostimulation system is indicated for SCS as an aid in the management of the following conditions:

- chronic, intractable pain of the trunk and/or limbs.
- stable intractable Angina Pectoris in patients who are not candidates for revascularization.*
- stable intractable Peripheral Vascular Disease of Fontaine Stage II or higher in patients who are not candidates for revascularization.*

* Physicians using SCS to treat these conditions are encouraged to share their results in scientific exchange or with Medtronic.

Intended use

The Medtronic SCS system delivers electrical stimulation to selected areas of the spinal cord to treat indicated patients. The Medtronic SCS system consists of a multi-programmable neurostimulator, one or more leads, and one or more optional extensions to connect the leads to the neurostimulator. An external clinician programmer is used to noninvasively adjust patient stimulation parameters. Patients may also make stimulation adjustments, based on their pain needs and clinician-prescribed programming options, via an external patient programmer.

Additional warning

In addition to the warnings summarized in the Information for Prescribers (IFP) booklet, the implantation of an SCS system for angina includes the following warning:

Intractable angina pectoris - Patients treated for intractable angina pectoris should be educated on the signs and symptoms of myocardial infarction and should seek medical attention immediately if signs and symptoms develop.

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
As a registered Notary Public for the State of Minnesota,
 I hereby certify that this is a true and accurate copy of the
 original document. The original is filed at Medtronic, Inc.

Indications sheet • Brochure concernant les indications • Indikationenblatt •
 Hoja de indicaciones • Indicatieboekje • Foglietto illustrativo •
 Vedlegg om indikationer • Käyttöaihelehtinen • Lista över indikationer •
 Φύλλο ενδείξεων • Indikationsside • Folheto de indicações •
 Folha de indicações • Informacje o wskazaniach • Informace o indikacích •
 Javallatokat ismertető lap • Перечень показаний • Zoznam indikácií •
 Endikasyonlar sayfası • Listă de indicații • Лист с показания •
 List sa indikacijama • Priloga z indikacijami • List s indikacijama •
 ورقة دواعي الاستعمال • גיליון התוויות

Ashley M. Luloff
 November 15, 2019

Indications

Neurostimulation for Peripheral Nerve Stimulation (PNS) using percutaneous leads - A Medtronic PNS neurostimulation system is indicated for PNS as an aid in the management of chronic, intractable pain of the posterior trunk.

Neurostimulation for Peripheral Nerve Stimulation (PNS) using surgical leads - A Medtronic PNS neurostimulation system is indicated for PNS as an aid in the management of chronic, intractable pain of the trunk and/or limbs.

Intended use

The Medtronic PNS system delivers electrical stimulation to selected areas of the peripheral nervous system to treat indicated patients. The Medtronic PNS system consists of a multi-programmable neurostimulator, one or more leads, and one or more optional extensions to connect the leads to the neurostimulator. An external clinician programmer is used to noninvasively adjust patient stimulation parameters. Patients may also make stimulation adjustments, based on their pain needs and clinician-prescribed programming options, via an external patient programmer.

Additional precautions for surgical leads

In addition to the precautions summarized in the Information for Prescribers (IFP) booklet, the implantation of a peripheral nerve stimulation system also includes the following precautions:

Peripheral nerve, handling - Be extremely careful to avoid excess manipulation and irritation of the peripheral nerve. If the nerve is irritated, the patient may experience pain, dysesthesia, or uncomfortable stimulation.

Peripheral nerve, location - Avoid implanting the paddle and the lead body in an area of repetitive movement (eg, over a joint). Repetitive movement can break a lead wire, resulting in intermittent or loss of stimulation.

Additional adverse event

In addition to the adverse events summarized in the Information for Prescribers (IFP) booklet, the implantation of a peripheral nerve stimulation system also includes the following risk:

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
MEDTRONIC CONTAINS A REPRODUCTION OF THE MEDTRONIC LOGO. IT IS NOT TO BE REPRODUCED OR USED IN ANY MANNER WITHOUT THE WRITTEN PERMISSION OF MEDTRONIC.


LASER

NOTE: IMAGES NOT TO SCALE

PART NUMBER
A93861
REV
A
SHEET
3 of 3


Specify™ SureScan™ MRI 5-6-5 977C1XX

REF 977C1XX
Reorder number

LOT WWW000000W
Lot number

Use by YYYY-MM-DD
Date of manufacture

EC REP
E.C. Representative
Medtronic B.V.
Earl Bakkerstraat 10
6422 PJ Heerlen
The Netherlands

www.medtronic.com

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Villalba, Puerto Rico, USA
Product of USA

Medtronic
977C1XX

LOT WWW000000W
Lot number

Use by YYYY-MM-DD
Date of manufacture

PIN W000WWW00W

Lead kit

XX cm
Length

MB04321A001 A - MB04325A001 A

SHIELD BOX LABEL



Medtronic
977C1XX

LOT WWW000000W
Lot number

Use by YYYY-MM-DD
Date of manufacture

Lead kit

Specify™ SureScan™ MRI 5-6-5 977C1XX


Length XX cm

REF 977C1XX
Reorder number

LOT WWW000000W
Lot number

Use by YYYY-MM-DD
Date of manufacture

Medtronic
977C1XX

LOT WWW000000W
Lot number

Use by YYYY-MM-DD
Date of manufacture

EC REP
E.C. Representative
Medtronic B.V.
Earl Bakkerstraat 10
6422 PJ Heerlen
The Netherlands

www.medtronic.com

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Villalba, Puerto Rico, USA
Product of USA

Medtronic
977C1XX

LOT WWW000000W
Lot number

Use by YYYY-MM-DD
Date of manufacture

PIN W000WWW00W

MB04321A001 A - MB04325A001 A

STERILE PACK LABEL

OPENING INSTRUCTION LABEL

(01) 00000000000000 (17) YYMMDD (21) WWW000000W
FOR YOUR RECORDS

Medtronic 977C1XX	LOT WWW000000W Lot number	Medtronic 977C1XX	LOT WWW000000W Lot number
Medtronic 977C1XX	LOT WWW000000W Lot number	Medtronic 977C1XX	LOT WWW000000W Lot number
Medtronic 977C1XX	LOT WWW000000W Lot number	Medtronic 977C1XX	LOT WWW000000W Lot number
Medtronic 977C1XX	LOT WWW000000W Lot number	Medtronic 977C1XX	LOT WWW000000W Lot number

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*Ashley M. Luloff
November 15, 2019*

ASHLEY M. LULOFF
NOTARY PUBLIC - MINNESOTA
MY COMMISSION EXPIRES 01/31/24

LASER

NOTE: IMAGES NOT TO SCALE

PART NUMBER
A93863
REV
B
SHEET
3 of 3


Specify™ SureScan™ MRI 2x8 977C2XX

Neurostimulation for Chronic Pain Lead kit MR Conditional

XX cm Length

2.8 mm

REF 977C2XX
Reorder number

LOT WWW000000W
Lot number

Use by
YYYY-MM-DD
Date of manufacture

Package contents Lead Product documentation Accessories

Do not sterilize Do not re-use Consult instructions for use

EC REP
E.C. Representative
Medtronic B.V.
Earl Bakkerstraat 10
5422 PJ Heerlen
The Netherlands

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in:
Villalba, Puerto Rico, USA
Product of USA

CE0123

Medtronic
977C2XX

PIN W000WWW00W

Lead kit

XX cm Length

LOT WWW000000W
Lot number

Use by
YYYY-MM-DD
Date of manufacture

M864321A01 A - M864327A01 B

Specify™ SureScan™ MRI 2x8 977C2XX

Neurostimulation for Chronic Pain Lead kit MR Conditional

2.8 mm

XX cm Length

REF 977C2XX
Reorder number

LOT WWW000000W
Lot number

Use by
YYYY-MM-DD
Date of manufacture

Package contents Lead Accessories

Do not sterilize Do not re-use

EC REP
E.C. Representative
Medtronic B.V.
Earl Bakkerstraat 10
5422 PJ Heerlen
The Netherlands

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in:
Villalba, Puerto Rico, USA
Product of USA

CE0123

Medtronic
977C2XX

PIN W000WWW00W


Lead kit

XX cm Length

LOT WWW000000W
Lot number

Use by
YYYY-MM-DD
Date of manufacture

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XX cm Length

LOT WWW000000W
Lot number

Use by
YYYY-MM-DD
Date of manufacture

M864321A01 A - M864327A01 B

977C2XX

Lead kit

XX cm Length

LOT WWW000000W
Lot number

Use by
YYYY-MM-DD
Date of manufacture

M864321A01 A - M864327A01 B

SHIELD BOX LABEL STERILE PACK LABEL

SIDE LABEL OPENING DISTRIBUTION LABEL

FOR YOUR RECORDS

Medtronic 977C2XX Lot number	Medtronic 977C2XX Lot number
Medtronic 977C2XX Lot number	Medtronic 977C2XX Lot number
Medtronic 977C2XX Lot number	Medtronic 977C2XX Lot number
Medtronic 977C2XX Lot number	Medtronic 977C2XX Lot number

ASHLEY M LULOFF
NOTARY PUBLIC - MINNESOTA
MY COMMISSION EXPIRES 01/31/24

PISCES-QUAD™ PLUS

3888-XX

Neurostimulation
for Chronic Pain

Lead kit

XX cm

Length


REF 3888-XX
Reorder number

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

2000-00-00
YYMMDD
Date of manufacture

34 °C
-30 °F

57 °C
135 °F

STERILE EO

Do not
re-use

Consult
instructions
for use

Sterilized using
ethylene oxide

EC REP
E.C. Representative
Medtronic B.V.
Earl Bakkenstraat 10
8422 PJ Heerlen
The Netherlands

www.medtronic.com
/manuals

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in:
Vila Real, Portugal, USA
Product of USA

CE0123

3888-XX

Lead kit

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

XX cm
Length

PIN W000WWW00W

M06147A001 A - M061562A001 A

SHELF BOX LABEL

BAR CODE LABEL

3888-XX

Lead kit

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

M06147A001 A - M061562A001 A

SIDE LABEL

FOR YOUR RECORDS

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

Medtronic
3888-XX
LOT WWW000000W
Lot number

PISCES-QUAD™ PLUS


3888-XX

Neurostimulation for
Chronic Pain

Lead kit

Length XX cm


STERILE

REF 3888-XX
Reorder number

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

2000-00-00
YYMMDD
Date of manufacture

34 °C
-30 °F

57 °C
135 °F

STERILE EO

Do not
re-use

Consult
instructions
for use

Sterilized using
ethylene oxide

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

Do not use
if package
is damaged

Consult
instructions
for use

STERILE EO

Sterilized using
ethylene oxide

www.medtronic.com
/manuals

EC REP
E.C. Representative
Medtronic B.V.
Earl Bakkenstraat 10
8422 PJ Heerlen
The Netherlands

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Vila Real, Portugal, USA
Product of USA

CE0123


PIN W000WWW00W

M06147A001 A - M061562A001 A

STERILE PACK LABEL

OPENING INSTRUCTION LABEL

* 220290001 *

NEURO LABEL SET
STIM STANDARD
220290001 (DIE B-10526001)

ASHLEY M LUILOFF
NOTARY PUBLIC - MINNESOTA
MY COMMISSION EXPIRES 01/31/24

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PISCES-QUAD™

3487A-XX

Neurostimulation for Chronic Pain

Lead kit

XX cm
Length

4.3mm

REF 3487A-XX
Reorder number

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

2000-00-00
YYMMDD
Date of manufacture

Package contents

Lead

Product documentation

Accessories

57 °C
135 °F

-34 °C
-30 °F

Temperature limitations

Do not re-use

Consult instructions for use

STERILE EO

Sterilized using ethylene oxide

EC REP

E.C. Representative
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands

www.medtronic.com /manuals

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Vila Rica, Puerto Rico, USA
Product of USA

CE0123

Medtronic

3487A-XX

Lead kit

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

XX cm
Length

PIN W000WW00W

M061475A001 A - M061551A001 A

SHIELD BOX LABEL

Medtronic

PISCES-QUAD™

3487A-XX

Neurostimulation for Chronic Pain

Lead kit

4.3mm

XX cm
Length

STERILE

REF 3487A-XX
Reorder number

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

2000-00-00
YYMMDD
Date of manufacture

Package contents

Lead

Accessories

57 °C
135 °F

-34 °C
-30 °F

Temperature limitations

Do not re-use

Consult instructions for use

STERILE EO

Sterilized using ethylene oxide

EC REP

E.C. Representative
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands

www.medtronic.com /manuals

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Vila Rica, Puerto Rico, USA
Product of USA

CE0123

Medtronic

3487A-XX

Lead kit

LOT WWW000000W
Lot number


1975-09-17
YYMMDD
Use by

XX cm
Length


PIN W000WW00W

M061475A001 A - M061551A001 A

STERILE PACK LABEL

NEURO LABEL SET
STIM STANDARD
220200001 (DIE B-105326001)

Ashley M. Luloff November 15, 2019.

Medtronic

3487A-XX

Lead kit

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

XX cm
Length

PIN W000WW00W

M061475A001 A - M061551A001 A

SHIELD BOX LABEL

Medtronic

3487A-XX

Lead kit

LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by

XX cm
Length

PIN W000WW00W

M061475A001 A - M061551A001 A


STERILE PACK LABEL

FOR YOUR RECORDS

Medtronic 3487A-XX Lot number LOT WWW000000W	Medtronic 3487A-XX Lot number LOT WWW000000W
Medtronic 3487A-XX Lot number LOT WWW000000W	Medtronic 3487A-XX Lot number LOT WWW000000W
Medtronic 3487A-XX Lot number LOT WWW000000W	Medtronic 3487A-XX Lot number LOT WWW000000W
Medtronic 3487A-XX Lot number LOT WWW000000W	Medtronic 3487A-XX Lot number LOT WWW000000W

As a registered Notary Public for the State of Minnesota, I hereby certify that this is a true and accurate copy of the original document. The original is filed at Medtronic, Inc.

PART NUMBER	REV	SHEET
A78004	B	2 of 2

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





NOTARY PUBLIC - MINNESOTA

MY COMMISSION EXPIRES 01/31/24

As a registered Notary Public for the State of Minnesota, I hereby certify that this is a true and accurate copy of the original document. The original is filed at Medtronic, Inc.

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
As a registered Notary Public for the State of Minnesota, I hereby certify that this is a true and accurate copy of the original document. The original is filed at Medtronic, Inc.



Vectris™ SureScan™ MRI 977A2XX

REF 977A2XX
Reorder number
LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by
2000-00-00
YYMMDD
Date of manufacture

EC REP
E.C. Representative
Medtronic B.V.
Eet Balkenstraat 10
6422 PJ Heerlen
The Netherlands

www.medtronic.com /manuals

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Vialba, Puerto Rico, USA
Product of USA



PIN W000WWW00W
M061475A001 A - M061383A001 A

SHELF BOX LABEL



Medtronic	LOT WWW000000W
977A2XX	Lot number
Medtronic	LOT WWW000000W
977A2XX	Lot number
Medtronic	LOT WWW000000W
977A2XX	Lot number
Medtronic	LOT WWW000000W
977A2XX	Lot number

Vectris™ SureScan™ MRI 977A2XX

REF 977A2XX
Reorder number
LOT WWW000000W
Lot number

1975-09-17
YYMMDD
Use by
2000-00-00
YYMMDD
Date of manufacture


EC REP
E.C. Representative
Medtronic B.V.
Eet Balkenstraat 10
6422 PJ Heerlen
The Netherlands

www.medtronic.com /manuals

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
www.medtronic.com

Manufactured in
Vialba, Puerto Rico, USA
Product of USA


PIN W000WWW00W
M061475A001 A - M061383A001 A

STERILE PACK LABEL

Medtronic	LOT WWW000000W
977A2XX	Lot number
Medtronic	LOT WWW000000W
977A2XX	Lot number
Medtronic	LOT WWW000000W
977A2XX	Lot number
Medtronic	LOT WWW000000W
977A2XX	Lot number

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Medtronic

MRI guidelines for Medtronic neurostimulation systems for chronic pain
See "START HERE" section before conducting MRI.

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
USA Rx only

CE0123

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic, Inc.

Ashley M. Luloff
November 15, 2019

Explanation of symbols on product or package labeling
Refer to the appropriate product for symbols that apply.

CE0123 Conformité Européenne (European Conformity). This symbol means that the device fully complies with AAMI Directive 90/385/EEC (NB 0123).

Manufacturer

Authorized Representative in the European Community

For USA audiences only

Magnetic Resonance (MR) Conditional

Magnetic Resonance (MR) Unsafe

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. Invi™ and SureScan™ are trademarks of a Medtronic company.

4 English NPI guidelines for Medtronic neurostimulation systems for chronic pain 2017-08-01

Table of contents

Introduction	7
Neurostimulator model numbers	7
Patient ID card	7
Obtain the latest MRI guidelines labeling	8
Clinician programmer and patient control device	8
General information on MRI procedures and neurostimulation system interactions	8
Types of electromagnetic fields generated by MRI systems	8
Potential interactions for implanted neurostimulation systems in the MRI environment	9
START HERE – Eligibility Identification	10
Use the identification checklist to identify the patient's scan eligibility first	10
Warnings	10
Precautions	10
Identification checklist	10
Full-body eligible MRI scan conditions	14
Full-body eligible – MRI equipment and scan requirements	14
Full-body eligible – Preparing the patient before the MRI scan	16
Full-body eligible – Pre-MRI scan operations and considerations	17
Image artifacts and distortion	17
Full-body eligible – During the MRI scan	18
Full-body eligible – Post-MRI scan	19
Head-only eligible MRI scan conditions	21
Checklist before proceeding with the head-only scan conditions	21
Head-only eligible – MRI equipment and scan requirements	23
Head-only eligible – Preparing the patient before the MRI scan	25
Head-only eligible – Pre-MRI scan operations and considerations	25
Head-only eligible – During the MRI scan	26
Head-only eligible – Post-MRI scan	26
Appendix A: Examples of MRI scan-type eligibility sheets	28
Examples of MRI reports and screen printouts	28
Appendix B: MRI Scan-Type Eligibility Form	31
Medtronic Neurostimulation System MRI Scan-Type Eligibility	31

5 English MRI guidelines for Medtronic neurostimulation systems for chronic pain 2017-06-01

Introduction

It is important to read the information in this manual in its entirety before conducting a magnetic resonance imaging (MRI) examination on a patient with any implanted component of a Medtronic neurostimulation system for chronic pain.

Contact a Medtronic representative if you have any questions about the information in this manual.

Neurostimulator model numbers

The neurostimulator model numbers listed herein are MRI Conditional.

Do not use model numbers alone to determine which MRI scan conditions to use in these MRI guidelines. Always begin with the "START HERE – Eligibility Identification" section of this manual and use the yes/no identification checklist that starts on page 10 to determine the patient's MRI scan-type eligibility and the appropriate scan conditions to use for the patient's implanted Medtronic neurostimulation system for chronic pain.

Follow these MRI guidelines and conditions for approved indications to determine whether and how to perform an MRI scan safely on a patient with a fully implanted Medtronic neurostimulation system for chronic pain with the neurostimulator model number listed herein.

These MRI guidelines apply to the following Medtronic implanted neurostimulator model numbers:

97702	37701	37711	7479
97712	37702	37712	74708
97713	37703	37713	7427
97714	37704	37714	7427V
97715			7426 (see the warning)
97716			

Warning: Medtronic recommends physicians not prescribe MRI for a patient who has an implanted level 3 Model 7426 Neurostimulator. The level 3 Neurostimulator is highly susceptible to reset or damage when subjected to an MRI examination. If reset, the neurostimulator must be reprogrammed. If damaged, the neurostimulator must be replaced. The level 3 Neurostimulator has an increased risk of induced electrical current, which may stimulate or shock the patient.

Patient ID card

Advise the patient to bring the most up-to-date patient ID card to all MRI appointments. MRI personnel can then use the patient ID card to identify Medtronic as the manufacturer.

of the patient's neurostimulation system and to confirm the model number of the implanted neurostimulator.

Obtain the latest MRI guidelines labeling

Always obtain the latest MRI guidelines. Refer to the contact information at the back of this manual, or go to www.medtronic.com/mri.

Copies of these MRI guidelines may not be the most up-to-date version if not received directly from the website or in another manner from Medtronic the same day of the patient's MRI appointment.

Clinician programmer and patient control device

For Medtronic neurostimulation systems with SureScan MRI Technology, external control devices (ie, a clinician programmer or a patient control device) are used to determine MRI scan-type eligibility and are used to place the neurostimulation system in MRI mode, which turns stimulation off. Inform the patient with a neurostimulation system that stimulation needs to be turned off prior to the MRI scan.

If the patient brought a patient control device to the MRI appointment, go to "START HERE – Eligibility Identification" on page 10 and use the identification checklist in that section.

If the clinician programmer or a patient control device cannot communicate with the implanted neurostimulation system or if the neurostimulator has reached EOS (end of service), then MRI scan-type eligibility cannot be confirmed via the external control devices. Resynchronizing the implanted neurostimulation system configuration from the patient's medical records is required. Unless the implanted system configuration is known and it is determined to be safe to perform an MRI under specific conditions, an MRI scan should not be conducted.

For operation of the clinician programmer, refer to the appropriate clinician programmer software manual for those instructions.

General information on MRI procedures and neurostimulation system interactions

Types of electromagnetic fields generated by MRI systems

An MRI system produces 3 types of electromagnetic fields that may interact with implanted device systems. All 3 of these fields are necessary to produce an MRI image. The 3 fields are defined as follows:

Static magnetic field – This is a steady state non-varying magnetic field that is always present around an MRI machine, even when no scan is underway.

Gradient magnetic fields – These low-frequency pulsed magnetic fields are present only during a scan. MRI equipment uses 3 orthogonal gradient magnetic fields to construct the 3-dimensional image.

RF field – This is a pulsed radio-frequency (RF) field that is present only during a scan. The RF field can be produced by a variety of transmission RF coils, such as a whole body transmit coil (that is built into the scanner) or an extremity coil (for example, a transmit/receive head coil).

Potential interactions for implanted neurostimulation systems in the MRI environment

The Medtronic neurostimulation systems with SureScan MRI Technology have been designed to minimize the potential interactions described in this section when the appropriate conditions described in this manual are followed.

Heating – The RF fields generated by an MRI scanner induce RF energy onto an implanted lead system that may cause heating at the lead electrodes or along the lead body. In addition, the gradient magnetic and RF fields may cause heating of the neurostimulator.

Note: Heating can occur even if only a lead or extension is implanted.

Factors that increase the risks of heating and patient injury include, but are not limited to, the following:

- High MRI specific absorption rate (SAR) RF power levels
- Low impedance leads or extensions (Medtronic product names or model numbers designated by a "Z," an "LZ," or "low impedance")
- Implanted lead systems with small surface area electrodes
- Short distances between lead electrodes and heat-sensitive tissue

Magnetic field interactions – The magnetic material of an implanted system may exert force, vibration, and torque effects due to the static magnetic field and gradient magnetic fields produced by an MRI scanner. Patients may feel a mild tugging or vibration sensation at the site of the device implant. Patients being scanned with recent implant incisions should be monitored for any surgical wound discomfort.

Induced stimulation – The gradient magnetic and RF fields produced by an MRI scanner induce energies onto an implanted lead system that could potentially cause unintended stimulation, which the patient could experience as a tingling, shocking, or jolting sensation.

Note: Induced stimulation can occur even if only a lead or extension is implanted.

Device damage – The voltages induced by the MRI fields may damage the neurostimulator electronics requiring reprogramming, explantation, or replacement.

Device interactions – MRI may affect the operation of the neurostimulator and require reprogramming of the neurostimulator with the clinician programmer after the MRI scan. The MRI may also reset the parameters to power-on-reset (POR) settings, which may also require reprogramming of the neurostimulator after the MRI scan.

START HERE – Eligibility identification

Use the identification checklist to identify the patient's scan eligibility first

Use the yes/no identification checklist that starts on page 10 in this section to determine the patient's MRI scan-type eligibility and the appropriate scan conditions to use for the patient's implanted Medtronic neurostimulation system for chronic pain.

The MRI scan-type eligibility depends on a combination of factors pertaining to the patient's implanted neurostimulation system.

Warnings

Other implanted devices – Prior to an MRI examination, determine whether the patient has multiple medical device implants, either active medical device implants (such as deep-brain stimulation systems, implantable cardiac defibrillators, etc) or passive medical device implants (such as spinal hardware, stents, etc). The most restrictive MRI exposure requirements must be used of the medical device implants. Contact the appropriate device manufacturers if you have questions. If you are unclear what implants may be present, perform an x-ray to determine implant type and location. Do not conduct an MRI examination if any conditions or implants that would prohibit or contraindicate an MRI are present.

Trial systems (neurostimulation systems that are not fully implanted) – Physicians should not prescribe MRI for patients undergoing trial stimulation or who have any neurostimulation system components that are not fully implanted. Explant all trial stimulation components if an MRI scan is required. MRI has not been tested on trial stimulation components and may cause heating of the lead electrodes, resulting in tissue damage or serious patient injury.

Precautions

External devices are MR Unsafe in the scanner (magnet) room – Do not allow the following Medtronic external control devices into the MRI scanner (magnet) room. These devices contain ferromagnetic material, which can be affected by the MRI magnet. These devices are MR Unsafe:

- Patient control device
- Recharger
- External neurostimulator
- Clinician programmer


Identification checklist

1. Did you receive a sheet denoting the MRI scan-type eligibility for the patient's MRI appointment?

See "Appendix A: Examples of MRI scan-type eligibility sheets" on page 28.

- ☐ **Yes** (1) Confirm the patient's name and date on the eligibility sheet. The date on the eligibility sheet should be on or near the date of the MRI appointment.
- Note:** The further the date on the eligibility sheet is from the patient's MRI appointment, the greater the chance that the following occurred:
- The patient had an event (eg, revision surgery of the implanted neurostimulation system) that may have changed the scan eligibility.
 - The patient's stimulation was turned back on.
- (2) Go to step 4 on page 12.
- ☐ **No** Go to the next step.
-

2.

Model number is on the back

Patient control device
(type A)

Model number is on the front

Patient control device
(type B)

Did the patient bring a Model 977__ series patient control device to the MRI appointment?

- ☐ **Yes** (1) Ask the patient to navigate to the MRI Mode screen to activate MRI mode.
- **MRI mode:** When MRI mode is activated:
 - The screen on the patient control device displays MRI scan-type eligibility.
 - If the eligibility screen displays, it means stimulation has been turned off.
 - If the patient is unable to navigate to the MRI Mode screen, contact Medtronic Technical Services for assistance.
 - Do not deactivate, or exit, MRI mode or turn stimulation on with the patient control device until after the patient's MRI scan is complete and the patient is outside of the scanner (magnet) room.
- (2) After the MRI Mode screen displays, go to step 4 on page 12.
- ☐ **No** Go to the next step.
-

3. Does the patient have a Medtronic patient ID card with a neurostimulator model number that matches one of the neurostimulator model numbers listed on page 7?

- ☐ Yes Go to "Head-only eligible MRI scan conditions" on page 21.
- ☐ No The patient's implanted system cannot be determined. Stop. These MRI guidelines do not apply. Go to www.medtronic.com/mri or call Medtronic Technical Services.

4.

MRI Conditional Full Body Scan Eligible

Does the MRI scan-type eligibility sheet or patient control device denote (with text and/or all of these symbols) full-body scan eligible?

- ☐ Yes (1) Confirm that the neurostimulator model number on the eligibility sheet or on the MRI Mode screen of the patient control device matches one of the model numbers on page 7.
- (2) If yes, go to "Full-body eligible MRI scan conditions" on page 14.
- ☐ No Go to the next step.

5.

MRI Conditional Head Scan Eligible with Transmit/Receive Head Coil

Does the MRI scan-type eligibility sheet or patient control device denote (with text and/or all of these symbols) head-only scan eligible?

- ☐ Yes (1) Confirm that the neurostimulator model number on the eligibility sheet or on the MRI Mode screen of the patient control device matches one of the model numbers on page 7.
- (2) Go to "Head-only eligible MRI scan conditions" on page 21.
- ☐ No Go to the next step.

6.

The neurostimulation system MRI scan-type eligibility cannot be determined.

Does the MRI scan-type eligibility sheet or patient control device denote (with text and/or all of these symbols) that the MRI scan-type eligibility cannot be determined?

- ☐ **Yes**
- (1) Confirm that the neurostimulator model number on the eligibility sheet or on the MRI Mode screen of the patient control device matches one of the model numbers on page 7.
 - (2) Go to "Head-only eligible MRI scan conditions" on page 21.
- ☐ **No** Stop. These MRI guidelines do not apply. Go to www.medtronic.com/mri or call Medtronic Technical Services.
-

Notes:

- The "consult instructions for use" symbol () when shown with MRI scan eligibility means "consult the MRI guidelines for this neurostimulation system."
- If a patient control device was used to determine eligibility, make a photocopy of the MRI Mode screen (shows the scan-type eligibility) that displayed on the patient control device, if needed.
- For interpretation of the information code () if given on the eligibility sheet or on the MRI Mode screen of the patient control device, call Medtronic Technical Services.
- Do not deactivate, or exit, MRI mode or turn stimulation on with the patient control device until after the patient's MRI scan is complete and the patient is outside of the scanner (magnet) room.

Full-body eligible MRI scan conditions

MR Conditional Full Body Scan Eligible

Before proceeding with this full-body eligible section, confirm that the "START HERE – Eligibility Identification" section (starts on page 10) has been followed and full-body scan eligibility has been determined by either an MRI scan-type eligibility sheet or a patient control device.

A patient who is "MRI-CS (or MR Conditional) full body scan eligible" can have any part of the anatomy scanned when the specific conditions in this full-body eligible section are met.

Full-body eligible – MRI equipment and scan requirements

Table 1. Full-body eligible – MRI equipment and scan requirements

MRI system type	1.5-T horizontal closed bore with maximum spatial gradient of 19 T/m (1900 gauss/cm). ⚠ Warning: Only use 1.5-T horizontal closed bore MRI systems. Other MRI systems (such as 0.6-T or 3.0-T, and open bore machines) have not been tested and could cause device damage and excessive heating, which can result in tissue damage or serious patient injury.
MRI manufacturers	No restrictions.
Radio-frequency (RF) frequency	Approximately 64 MHz. ⚠ Warning: Do not conduct MRI scans with nonproton scanning frequencies (such as, ¹³ C, ²³ Na, or ³¹ P). Frequencies other than 64 MHz have not been tested and could cause device damage and excessive heating, which can result in tissue damage or serious patient injury.
RF coils	Transmit coil: - quadrature transmit/receive body coil (built-in). - quadrature transmit/receive head coil. ⚠ Warning: Do not use RF transmit coils other than a body transmit/receive (built-in) quadrature coil or a head transmit/receive quadrature coil. Other transmit/receive coils (eg, linear coils) have not been tested and could cause excessive heating, which can result in tissue damage or serious patient injury. Receive-only coil: any type.

Table 1. Full-body eligible – MRI equipment and scan requirements (continued)

RF power	Use Normal Operating Mode. Specific absorption rate (SAR): - Whole body SAR: must be ≤ 2.0 W/kg (0.91 W/lb) as reported by the MRI equipment. - Head SAR: must be ≤ 3.2 W/kg (1.45 W/lb) as reported by the MRI equipment.  Warning: Do not conduct MRI scans in the following modes: - First Level Controlled Operating Mode - Second Level Controlled Operating Mode (i.e. research mode) These modes allow higher levels of RF energy and may cause excessive heating, which can result in tissue damage or serious patient injury.
Gradients	Gradient systems with a maximum gradient slew rate performance per axis of 200 T/m/s or less.  Warning: Do not use gradient systems producing gradient slew rates greater than 200 T/m/s because they have not been tested and could cause increased risk of induced stimulation (resulting in shocking or jolting sensations, discomfort, or pain for the patient) or heating of the neurostimulator.
Active scan time limits	MRI scan durations should not exceed a total of 30 minutes of active scan time within a 90-minute window (within every 90-minute window should be a total of 60 minutes of nonscan time).  Warning: Do not exceed a total of 30 minutes of active scan time within a 90-minute window. Exceeding the active scan time duration increases the risk of tissue heating.
Landmark (isocenter location)	No restrictions. All anatomical locations can be scanned.

Full-body eligible – Preparing the patient before the MRI scan

Table 2. Full-body eligible – Preparing the patient before the MRI scan

MRI mode on, stimulation off	 MRI Conditional Full Body Scan Eligible
	Placing the device in MRI mode turns stimulation off. The text and/or all of the symbols above denote full-body MRI scan eligibility and indicate that the implanted system is in MRI mode. Caution: Before conducting the MRI scan, confirm that the patient's implanted neurostimulation system is off. Leaving stimulation on during the scan could increase the potential for uncomfortable, unintended stimulation. If you are not certain if stimulation is off, ask the patient if it is off. Note: A depleted neurostimulator is considered off.
Core body temperature	Fever Warning: Do not perform an MRI scan if the patient's body temperature is above 38 °C (100 °F). Elevated body temperature in conjunction with tissue heating caused by an MRI scan increases the risk of excessive tissue heating, which could cause tissue damage. No blankets Warning: Do not cover the patient with blankets or heated blankets. Blankets raise the patient's body temperature and increase the risk of tissue heating, which could cause tissue damage.
Patient weight, minimum	Caution: Do not conduct MRI scans on patients who weigh less than 40 kg (88 lbs). Scans for patients who weigh less than 40 kg (88 lbs) have not been tested and could cause excessive heating, resulting in tissue damage.
Sedation	If possible, do not sedate the patient so that the patient can provide feedback of any problems during the examination.

Table 2. Full-body eligible – Preparing the patient before the MRI scan (continued)

Patient position within the bore	Position the patient in a prone or supine position in the MRI bore. ⚠ Warning: Do not position the patient in other positions, eg, on his or her side (called the lateral decubitus position) within the MRI bore. Scanning patients in positions other than prone or supine is untested and could cause excessive tissue heating during an MRI scan.
Inform the patient of risks	Inform the patient of all the risks of undergoing an MRI examination as stated in this full-body eligible section.
Patient communication with operator during scan	Instruct the patient to immediately inform the MRI operator if any discomfort, unexpected stimulation, shocking, or heating occurs during the examination.

Full-body eligible – Pre-MRI scan operations and considerations

Table 3. Full-body eligible – Pre-MRI scan operations and considerations

Enter patient weight	Enter the correct patient weight into the MRI console to ensure that the SAR is estimated correctly. ⚠ Warning: Ensure the patient weight is entered correctly to avoid the risk that the MRI scan is performed at an RF power level too high for the patient. An inappropriately high RF power level may cause excessive heating, which can result in tissue damage or serious patient injury.
Verify all parameters	Verify that all proposed MRI examination parameters comply with the MRI exposure requirements in this full-body eligible section. If not, the parameters must be modified to meet these requirements. If the parameters cannot be modified, do not perform an MRI.

Image artifacts and distortion

SureScan leads have demonstrated minimal image distortion for areas surrounding the implanted leads when the device is out of the field of view. Significant image distortion can result from the presence of the device within the field of view. Image artifacts and distortion resulting from the presence of the device and the leads within the field of view must be considered when selecting the field of view and imaging parameters. These factors must also be considered when interpreting the MRI images.

Careful choice of pulse sequence parameters, location of the angle, and location of the imaging plane may minimize MRI image artifacts. However, the reduction in image distortion obtained by adjustment of pulse sequence parameters will usually compromise signal-to-noise ratio.

The following general principles should be followed:

- Avoid using the body receive coil if possible. Use a local receive-only coil instead.
- Use imaging sequences with stronger gradients for both slice and read encoding directions. Use higher bandwidth for both radio-frequency pulse and data sampling.
- Choose an orientation for the read-out axis that minimizes the appearance of in-plane distortion.
- Use spin echo or gradient echo MR imaging sequences with a relatively high data sampling bandwidth.
- Use a shorter echo time for gradient echo technique, whenever possible.
- Be aware that the actual imaging slice shape can be curved in space due to the presence of the field disturbance of the neurostimulator.
- Identify the location of the implant in the patient, and when possible, orient all imaging slices away from the implanted neurostimulator.

⚠ Warnings:

- If the MRI targeted image area is near the neurostimulator, it may be necessary to move the neurostimulator to obtain an image, or use alternate imaging techniques. MRI images may be severely distorted or image target areas can be completely blocked from view near the implanted neurostimulation system components, especially near the neurostimulator.
- If the neurostimulator is removed, remove the entire neurostimulation system. Do not remove the neurostimulator and leave the lead system implanted as this can result in higher than expected lead heating. Excessive heating can result in tissue damage or serious patient injury.

Contact Medtronic Technical Services for more information about the expected extent and appearance of MRI image artifacts and distortion for various scan conditions.

Full-body eligible – During the MRI scan

Table 4. Full-body eligible – During the MRI scan

Keep track of active scan time	Keep track that the active scan time is within a 90-minute window. See "Active scan time limits" in Table 1.
Monitor the patient	Monitor the patient both visually and audibly. Check the patient between each imaging sequence. Discontinue the MRI examination immediately if the patient is unable to respond to questions or reports any problems.
Patient comfort	Heating may be felt at the neurostimulator site during the MRI scan. If the heating causes the patient discomfort, stop the MRI scan immediately. Consider applying an ice pack or cold compress to the location after the scan is stopped.

Table 4. Full-body eligible – During the MRI scan (continued)

Neurostimulator tugging, vibration	During the MRI scan, the patient may feel tugging and/or vibration of the neurostimulator. If the tugging or vibration causes the patient considerable discomfort, stop the MRI scan. If the neurostimulator is close to the MRI bore wall, consider using a pillow to keep the neurostimulator away from the bore wall to minimize vibration.
------------------------------------	---

Full-body eligible – Post-MRI scan

Table 5. Full-body eligible – Post-MRI scan

Patient feedback	Verify that the patient has not experienced adverse effects as a result of the MRI. Contact Medtronic to report any adverse effects.
------------------	--

Table 5. Full-body eligible – Post-MRI scan (continued)

Turn stimulation back on	After the scan has been completed, instruct the patient to see the clinician managing the patient's neurostimulation system to have the stimulation turned back on. Or, if the patient has brought a patient control device to the MRI appointment, instruct the patient (outside of the scanner room) to turn the stimulation back on using the patient control device. Notes: • The patient control device may need to be re-enabled to make it operable for the patient to turn stimulation back on.  For patient control device, type A: Tell the patient to do the following: 1. Press the Increase/Decrease key to wake up the patient control device. 2. Press and hold the Lock button on the Unlock screen. 3. Press the Exit MRI Mode button. Then the patient can turn on stimulation.  For patient control device, type B: Tell the patient to press the Sync  key. Then the patient can turn on stimulation. • Turning stimulation on takes the neurostimulator out of MRI mode. • If the patient control device cannot synchronize with the neurostimulator, or cannot turn stimulation back on, or displays a screen with the letters "POR" on it, instruct the patient to see the clinician managing the patient's neurostimulation system. Contact Medtronic to report the POR event.
---------------------------------	---

Head-only eligible MRI scan conditions

Before proceeding with this head-only eligible section, confirm that the "START HERE – Eligibility Identification" section (starts on page 10) has been followed, then proceed to the following checklist.

Checklist before proceeding with the head-only scan conditions

1.

MR Conditional Head Scan Eligible with Transmit/Receive Head Coil

Did the "START HERE – Eligibility Identification" section direct you to this head-only eligible section because the MRI scan-type eligibility sheet or screen on the patient control device denotes (with text and/or all of these symbols) head-only scan eligible?

- ☐ Yes Go to "Head-only eligible – MRI equipment and scan requirements" on page 23 and continue from there.
- ☐ No Go to the next step.

2. Confirm that no part of the implanted system (ie, neurostimulator, extensions, leads, abandoned leads) is within an RF transmit/receive head coil.

This can be confirmed with x-ray imaging of the neck and head region or referring to the patient's records, for example.

- ☐ Yes Yes, confirmed. Go to the next step.
- ☐ No No, could not confirm. Stop. Contact Medtronic Technical Services.

3.

The neurostimulation system MRI scan-type eligibility cannot be determined.

Did the "START HERE – Eligibility Identification" section direct you to this head-only eligible section because the MRI scan-type eligibility sheet or screen on the patient control device denotes (with text and/or all of these symbols) that the MRI scan-type eligibility cannot be determined?

- ☐ Yes Go to "Head-only eligible – MRI equipment and scan requirements" on page 23 and continue from there.

☐ No Go to the next step.

4. The patient brought only the patient ID card to the MRI appointment, and the neurostimulator model number on the card is listed on page 7.

Confirm that the neurostimulator stimulation parameters are set as follows:

Stimulation: Off
If you are not certain if stimulation is off, ask the patient if it is off.

Other parameters: No change

In addition for the Itriel 3 Model 7425 only:

Magnetic (reed) switch: Disabled

☐ Yes Yes, confirmed. Go to "Head-only eligible – MRI equipment and scan requirements" on page 23 and continue from there.

☐ No No, could not confirm. Contact Medtronic Technical Services.

If all of the instructions stated in this head-only eligible section are followed, MRI scans of the head only using an RF transmit/receive head coil may be safely performed.

Head-only eligible – MRI equipment and scan requirements

Table 6. Head-only eligible – MRI equipment and scan requirements

RF coils	Transmit/receive head coil only. Important: - The RF transmit-receive head coil must not cover any implanted system component. Implanted system components may be located per approved labeling and may be as close as 0 cm to the lower (caudad) edge of the head coil, but no part shall be inside the head coil. - Ensure that the RF transmit body coil is not used. - If you are unsure if your MRI system has RF transmit/receive head coil capability, consult the MRI manufacturer. Warnings: - An MRI examination of the head only (no other part of the body) can be conducted safely using an RF transmit/receive head coil when all instructions in this head-only eligible section are followed. - Do not place any part of the RF transmit/receive head coil over any implanted neurostimulation system component. If the head coil extends over any part of the patient's neurostimulation system, higher than normal heating may occur at the location of the implanted lead electrodes. In addition, if the patient's neurostimulation system has a broken lead wire and the head coil extends over any part of the patient's neurostimulation system, higher than normal heating may occur at the break or lead electrodes. Excessive heating can cause tissue damage or serious patient injury.
MRI system type	1.5-T horizontal closed bore with maximum spatial gradient of 19 T/m (1900 gauss/cm). Warning: Only use 1.5-T horizontal closed bore MRI systems. Other MRI systems (such as 0.6-T or 3.0-T, and open bore machines) have not been tested and could cause device damage and excessive heating, which can result in tissue damage or serious patient injury.
MRI manufacturers	No restrictions.

Table 6. Head-only eligible – MRI equipment and scan requirements (continued)

Radio-frequency (RF) frequency	Approximately 64 MHz. ⚠ Warning: Do not conduct MRI scans with nonproton scanning frequencies (such as, ^{13}C, ^{23}Na, or ^{31}P). Frequencies other than 64 MHz have not been tested and could cause device damage and excessive heating, which can result in tissue damage or serious patient injury.
RF power	Use Normal Operating Mode. Specific absorption rate (SAR): - Head SAR: must be $\leq 3.2 \text{ W/kg}$ (1.45 W/kg) as reported by the MRI equipment. ⚠ Warning: Do not conduct MRI scans in the following modes: - First Level Controlled Operating Mode - Second Level Controlled Operating Mode (ie, research mode) These modes allow higher levels of RF energy and may cause excessive heating, which can result in tissue damage or serious patient injury.
Gradients	Gradient systems with a maximum gradient slew rate performance per axis of 200 T/m/s or less. ⚠ Warning: Do not use gradient systems producing gradient slew rates greater than 200 T/m/s because they have not been tested and could cause increased risk of induced stimulation (resulting in shocking or jolting sensations, discomfort, or pain for the patient) or heating of the neurostimulator.
Active scan time limits	No restrictions.
Landmark (isocenter location)	Head only. Do not place any part of the RF transmit/receive head coil over any implanted neurostimulation system component.

Head-only eligible – Preparing the patient before the MRI scan

Table 7. Head-only eligible – Preparing the patient before the MRI scan

Confirm that stimulation is turned off	Confirm that the checklist on page 21 (at the beginning of this head-only eligible section) was followed. Caution: Before conducting the MRI scan, confirm that the patient's implanted neurostimulation system is off. Leaving stimulation on during the scan could increase the potential for uncomfortable, unintended stimulation. If you are not certain if stimulation is off, ask the patient if it is off. Note: A depleted neurostimulator is considered off.
Core body temperature	Fever No restrictions. Shivering No restrictions.
Patient weight, minimum	No restrictions.
Sedation	If possible, do not sedate the patient so that the patient can provide feedback of any problems during the examination.
Inform the patient of risks	Inform the patient of all the risks of undergoing an MRI examination as stated in this head-only eligible section.
Patient communication with operator during scan	Instruct the patient to immediately inform the MRI operator if any discomfort, unexpected stimulation, shocking, or heating occurs during the examination.

Head-only eligible – Pre-MRI scan operations and considerations

Table 8. Head-only eligible – Pre-MRI scan operations and considerations

Enter patient weight	Enter the correct patient weight into the MRI console to ensure that the head SAR is estimated correctly.
Verify all parameters	Verify that all proposed MRI examination parameters comply with the MRI exposure requirements in this head-only eligible section. If not, the parameters must be modified to meet these requirements. If the parameters cannot be modified, do not perform an MRI.

Head-only eligible – During the MRI scan

Table 9. Head-only eligible – During the MRI scan

Monitor the patient	Monitor the patient both visually and audibly. Check the patient between each imaging sequence. Discontinue the MRI examination immediately if the patient is unable to respond to questions or reports any problems.
Neurostimulator tugging, vibration	During the MRI scan, the patient may feel tugging and/or vibration of the neurostimulator. If the tugging or vibration causes the patient considerable discomfort, stop the MRI scan.

Head-only eligible – Post-MRI scan

Table 10. Head-only eligible – Post-MRI scan

Patient feedback	Verify that the patient has not experienced adverse effects as a result of the MRI. Contact Medtronic to report any adverse effects.
------------------	--

Table 18. Head-only eligible – Post-MRI scan (continued)

Turn stimulation back on After the scan has been completed, instruct the patient to see the clinician managing the patient's neurostimulation system to have the stimulation turned back on.
Or, if the patient has brought a patient control device to the MRI appointment, instruct the patient (outside of the scanner room) to turn the stimulation back on using the patient control device.

Notes:

- The patient control device may need to be re-enabled to make it operable for the patient to turn stimulation back on.

For patient control device, type A:

Tell the patient to do the following:

1. Press the Increase/Decrease key to wake up the patient control device.
2. Press and hold the Lock button on the Unlock screen.
3. Press the Exit MRI Mode button.

Then the patient can turn on stimulation.

For patient control device, type B:

Tell the patient to press the Sync key.

Then the patient can turn on stimulation.

- Turning stimulation back on takes the neurostimulator out of MRI mode for those neurostimulators with that feature.
- If the patient control device cannot synchronize with the neurostimulator, or cannot turn stimulation back on, or displays a screen with the letters "POR" on it, instruct the patient to see the clinician managing the patient's neurostimulation system. Contact Medtronic to report the POR event.

Appendix A: Examples of MRI scan-type eligibility sheets

This appendix shows examples of MRI scan-type eligibility sheets produced by Medtronic clinician programmers that are appropriate for the patient's MRI appointment.

The MRI scan-type eligibility sheet shows the scan-type eligibility for the patient's implanted neurostimulation system for chronic pain after the patient has seen the clinician managing the patient's neurostimulation system.

The following are the three types of MRI scan-type eligibility sheets:

- MRI Report
- MRI-CS Mode screen printout
- MRI Scan-Type Eligibility Form completed by the patient's physician (see "Appendix B: MRI Scan-Type Eligibility Form" on page 31)

Examples of MRI reports and screen printouts

See the following two pages for examples of MRI scan-type eligibility sheets produced by Medtronic clinician programmers that are appropriate for the patient's MRI appointment.

Medtronic neuromodulation
Restoreur THAMRI Model# 97712 ***
Patient John Doe
Session date 01/06/2010
Programmer 8840 SN WHF000238 107-07

MRI-CS Report

As of 01.06.2010 the neurostimulation system is:

HRI-CS FULL BODY SCAN ELIGIBLE

See labeling for safety conditions:
www.medtronic.com/srt

Confirm that stimulation is OFF prior to an MRI scan

Neurostimulator model number

07/206/7510
10:15
Revised: 11/2001 TRAMBI
MRI - CS Mode

The neuromuscular system
MRI - CS system cannot be
determined.

See labeling for safety
conditions:
www.medtronic.com/mri
00000000000
① Confirm that stimulation is OFF
prior to an MRI scan.

Not eligible for full body scan:

- ① Device Location
- ① Lead Model(s)
- ① Lead Tip Location(s)
- ① Extension Model(s)
- ① Abandoned Pain Lead(s)

MRSA 67712

2017-08-01 MFI guidelines for Medtronic neurostimulation systems for chronic pain English 29

Neurostimulator model number

Date

MRI Scan-Type Eligibility Report

Tony Martin 123-45-789

7775 ME*****

Medtronic NeuroModulation

Session Date:

Dec 22, 2015 3:14 PM

As of Dec 22, 2015

MR Conditional Head Scan Eligible with Transport/Remove Head Coil

See labeling for MRI scan conditions:
www.medtronic.com/mri

Confirm that stimulation is OFF prior to an MRI scan.

Not eligible for full body scan due to:
Lead Motion
Extensions Present

Information Code:1234567890

END OF REPORT

Dec 22, 2015 4:15 PM from Programmer 002-4567890

Example MRI Scan-Type Eligibility Report

Appendix B: MRI Scan-Type Eligibility Form

Medtronic Neurostimulation System MRI Scan-Type Eligibility

At the time of the MRI appointment:

1. See labeling for MRI scan conditions: www.medtronic.com/mri.
2. Confirm that stimulation is off prior to the MRI scan.

Patient name:	
Physician name, office, address, and phone number:	

Important: Use the MRI (or MRI-CS) mode on the clinician programmer or patient control device to enter the information and scan eligibility result below.		
Date and time eligibility was determined:	Neurostimulator model number:	Neurostimulator serial number:
☐		MR Conditional Full Body Scan Eligible
☐		MR Conditional Head Scan Eligible with Transmit/Receive Head Coil
☐		The neurostimulation system MRI scan-type eligibility cannot be determined.
① Information code:		(not applicable for FULL BODY SCAN ELIGIBLE)

32 English MPR guidelines for Medtronic neurostimulation systems for chronic pain 2017-06-01

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Numer modelu neurostymulatora

Medtronic NeuroModulation
RestoreULTRAMRI Model 9771 NMB*****
Patient John Doe
Session date 01/06/2010
Programmer 8840 SN NMF0000238-05-01
MRI-CS Report
As of 01/06/2010 the neurostimulation system is:
MRI-CS FULL BODY SCAN ELIGIBLE
See labeling for safety conditions:
www.medtronic.com/mri
Confirm that stimulation is OFF prior to an MRI scan

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Data

14:15
RestoreULTRAMRI
MRI-CS Mode

The neurostimulation system MRI-CS eligibility cannot be determined.

See labeling for safety conditions:
www.medtronic.com/mri
0000000000
Confirm that stimulation is OFF prior to an MRI scan.

Not eligible for full body scan:

- Device Location
- Lead Model(s)
- Lead Tip Location(s)
- Extension Model(s)
- Abandoned Pain Lead(s)

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Numer modelu neurostimulatora

Data

MRI Scan-Type Eligibility Report
Tony Martin 123-45-789
07715 ME*****

Medtronic Neurostimulation
Session Date:
Dec 22, 2015 2:14 PM

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Not eligible for full body scan due to:
Lead Model
Stimulation Present

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MRI-CS Report

As of 01/06/2010 the neurostimulation system is:

MRI-CS FULL BODY SCAN ELIGIBLE

See labeling for safety conditions:
www.medtronic.com/ori

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dátum

01/06/2010 14:15
RestoreULTRANE
MRI-CS Mode

**For neurostimulation system
MRI-CS eligibility cannot be
determined.**

**See labeling for safety
conditions**
www.medtronic.com/ori
**Confirm that stimulation is OFF
prior to an MRI scan.**

Not eligible for full body scan

- Device Location
- Lead Model(s)
- Lead Tip Location(s)
- Extension Model(s)
- Abandoned Pole Location

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dátum

MRI Scan-Type Eligibility Report
Tony Morin 123-45-789
77715-IME*****

Medtronic Neurostimulation
Session Date
Dec 22, 2015 3:14 PM

As of Dec 22, 2015

MRI Conditional Head Scan Eligible with Transcranial Receive Head Coil

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Not eligible for full body scan due to:
Lead Magnet
Stimulation Present

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Medtronic Neuromodulation
 Patient Model: 62712
 Patient: John Doe
 Session date: 08/08/2018
 Programmer: 0040 IN 0470000000-100-0
 Medtronic CS System
 As of 01:00:3040 the neuromodulation system is:
 MRI-CS FULL ON: SCAN 0470000000
 See labeling for safety conditions:
 Medtronic CS
 Confirm that stimulation is OFF prior to an MRI scan

Not eligible for full body scan:

- Device Location
- Lead location
- Lead tip location(s)
- External Magnet(s)
- Abandoned Pile Location

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MR Scan Type: MRS
Scan Name: 123-45-789
07/15/2015

Medtronic Neuromodulation
Session Date:
Dec 22, 2015 14 PM

As of Dec 22, 2015

MR Conditional Must Scan Signify with T weighted Resonance Head Coil

See listing for MR scan conditions:
www.medtronic.com/mr

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Not eligible for full body scan due to:
Lead/Module
Electronic Present

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Medtronic Neurostimulation
Restored TMRs: John Doe 07712 *****
Patient: John Doe
Session date: 01/06/2010
Programmer: 8840 SN NHF0000232 07712
MRI-CS Report
As of 01/06/2010 the neurostimulation
system is:
MRI-CS FULL BODY SCAN ELIGIBLE
See labeling for safety conditions:
www.medtronic.com/mri
***Confirm that stimulation is OFF
prior to an MRI scan***

01/06/2010

Celo modelu neurčitosti

WARNING                  

DO NOT                  

DO NOT                  

DO NOT                  

DO NOT                  

DO NOT                  

DO NOT                  

DO NOT               

we are right.

Číslo modelu neurostimulátora

Dátum

MR Scan-Type Safety Report
Form Number 123-05-100

Medtronic Neurostimulation
Scan Date
Dec 22, 2015 4:14 PM

As of Dec 22, 2015

MR Conditional Head Scan Safety with Transmit/Receive Head Coil

See labeling for MR scan conditions
www.medtronic.com/mr

Ensure that stimulation is OFF prior to an MR scan.

Not eligible for full body scan due to
Lead Material
Stimulation Present

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Nörostimölör model numarası

Medtronic neuromodulation
RestoreUL FRAME Model 9771 MB*****
Patient John Doe
Session date 01/06/2010
Programmer 8840 SN KHR000028-03-07
MRI-CS Report
As of 01/06/2010 the neurostimulation system is:
MRI-CS FULL BODY SCAN ELIGIBLE
See labeling for safety conditions:
www.Medtronic.com/mri
Confirm that stimulation is OFF prior to an MRI scan

Örnek MRI-CS raporu
(yukarıda metin)

Nörostimölör model numarası

Tarih

01/06/2010 16:15
RestoreUL FRAME
MRI-CS Mode
The neurostimulation system MRI-CS eligibility cannot be determined.
See labeling for safety conditions:
www.Medtronic.com/mri
0000000000
Confirm that stimulation is OFF prior to an MRI scan.
Not eligible for full body scan:
Device Location
Lead Mode(s)
Lead Tip Location(s)
Extension Model(s)
Abandoned Pain Lead(s)
MRI-CS Mode ekranı görüntüsü

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Nöroestimülatör model numarası

Tarih

MR Scan Type Eligibility Report
Form Number: 123-45-789
7778

Medronic Neuromodulation
Scan Date:
Dec 22, 2019 14:14 PM

As of Dec 22, 2019

MR Conditional Head Scan Eligible with Transmittal Receive Head Coil

See labeling for MR scan conditions
www.medtronic.com/mri

Confirm that stimulation is OFF prior to an MR scan.

Not eligible for full body scan due to
Lead Models
Stimulators Present

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Medtronic Neuromodulation
Restored TRAN1 Mode# 0712 Aug*****
Patient John Doe
Session date 01/06/2010
Programmer 8840 SN HKF0002136 03/07
NRI-CS Report
As of 01/06/2010 the neurostimulation
system is:
NRI-CS FULL BODY SCAN ELIGIBLE
See labeling for safety conditions:
www.Medtronic.com/mri
**Confirm that stimulation is OFF
prior to An MRI scan**

Data

Data

MTCN/PIN
18:15

Resistor/LI FRAME
MRI-3S Mode

The neurostimulation system
MRI-3S eligibility cannot be
determined.

See labeling for safety
conditions:
www.medtronic.com/mri
0000000000

Confirm that stimulation is OFF
prior to an MRI scan.

Not eligible for full body scan:

- ☐ Device Location
- ☐ Lead Model(s)
- ☐ Lead Tip Location(s)
- ☐ Extension Mode(s)
- ☐ Abandoned Pain Lead(s)

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Date

MR Scan-Type Eligibility Report Tony Martin 123-45-789 77715-123456789	Medtronic Neurostimulation Session Date: Dec 22, 2015 3:14 PM
  	
As of Dec 22, 2015 MR Conditional Head Scan Eligible with Transcranial Magnetic Head Coil See listing for MR scan conditions www.medtronic.com/mri Confirm that stimulation is OFF prior to an MRI scan. Not eligible for full body scan due to Lead Model Extension Present Information Code: 1234567890	
END OF REPORT Dec 22, 2015 4:15 PM from Programmer #P2000101234	

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Table of contents

Contraindications	5
Warnings	5
Precautions	5
Use in specific populations	8
Physician training	8
Storage and sterilization	8
System implant	8
Clinician programming	11
Patient programming and patient control devices	12
Charging system (for rechargeable neurostimulators)	12
Patient activities	13
Component disposal	14
Individualization of Treatment	14
Adverse events summary	15
Patient Counseling Information	16

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4 English Information for prescribers 2018-07-31

Contraindications

Diathermy - Do not use shortwave diathermy, microwave diathermy or therapeutic ultrasound diathermy (all now referred to as diathermy) on patients implanted with a neurostimulation system. Energy from diathermy can be transferred through the implanted system and can cause tissue damage at the location of the implanted electrodes, resulting in severe injury or death. Refer to Appendix A: Electromagnetic interference for further information.

Warnings

Electromagnetic interference (EMI) - Electromagnetic interference is a field of energy generated by equipment found in the home, work, medical or public environments that is strong enough to interfere with neurostimulator function. Neurostimulators include features that provide protection from electromagnetic interference. Most electrical devices and magnets encountered in a normal day are unlikely to affect the operation of a neurostimulator. However, sources of strong electromagnetic interference can result in the following:

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Refer to Table 1, starting on page 8, and Appendix A: Electromagnetic interference for information on sources of EMI, the effect of EMI on the patient and the neurostimulation system, and instructions on how to reduce the risk from EMI.

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- Allergic or immune system response to the implanted materials
- Infection
- Lead, extension, or neurostimulator erosion through the skin or migration
- Leakage of cerebrospinal fluid
- Loss of pain relief may return patients to their underlying pain condition
- Patients on anticoagulation therapies may be at greater risk for postoperative complications such as hematomas that can result in paralysis
- Persistent pain at the neurostimulator site
- Placement of the epidural lead-extension is a surgical procedure that may expose patients to risks of epidural hemorrhage, hematoma, or paralysis
- Radicular chest wall stimulation
- Seroma or hematoma at the neurostimulator site
- Change in stimulation, possibly related to cellular changes around the electrode(s), shifts in electrode position, loose electrical connections, lead or extension fractures, which has been described by some patients as uncomfortable stimulation (jolting or shocking sensation).
- Formation of reactive tissue around the lead in the epidural space can result in delayed spinal cord compression and paralysis, requiring surgical intervention. Time to onset can range from weeks to many years after implant.
- Stimulation-dependent gastrointestinal symptoms such as nausea, diarrhea, incontinence, or constipation
- Stimulation-dependent bladder symptoms such as urinary retention, incontinence, or frequency

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
STATE OF MINNESOTA


ANOKA COUNTY

As a responsible member of Medtronic Inc., I hereby attest that I believe the Дополнение к Эксплуатационной документации attached herein to be a true and accurate copy and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
Instruction for use

For

Pisces- Quad, Pisces- Quad Compact and Pisces- Quad Plus Leads

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On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Pisces- Quad, Pisces- Quad Compact and Pisces- Quad Plus Leads" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
Instruction for use

For


Vectris SureScan MRI 1x8 Subcompact and 1x8 Compact Leads

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On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Vectris SureScan MRI 1x8 Subcompact and 1x8 Compact Leads" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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2013

M939859A006 Rev A

2013-05
Printing instructions: Refer to the "Implant Manual"
category Table 1 in doc# A00002 for Neuro Core
European Printing Instructions.

M939859A006 Rev A

2013-05
Printing instructions: Refer to the "Implant Manual"
category Table 1 in doc# A00002 for Neuro Core
European Printing Instructions.

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
Instruction for use

For

Vectris 1x 8 Compact and 1x8 Subcompact Trial Screening Lead Kit

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On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Vectris 1x 8 Compact and 1x8 Subcompact Trial Screening Lead Kit" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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LOT

Lot number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with MDD 93/42/EEC (NB 0123).

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Preparing for surgery 11

Placing a percutaneous lead 11

Testing stimulation intraoperatively 13

Removing the needle and closing 13

Preparing components for trial screening 14

Removing components after trial screening 14

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6 English 977D160, 977D260 2012-11

Description

The Medtronic Vectris 1x8 Subcompact Model 977D160 and Vectris 1x8 Compact Model 977D260 Trial Screening Leads are part of a neurostimulation system for pain therapy.

Additional contraindication

The Medtronic Vectris 1x8 Subcompact Model 977D160 and Vectris 1x8 Compact Model 977D260 Trial Screening Leads are contraindicated for long-term implantation. These leads **MUST BE REMOVED** within thirty (30) days of implant.

Package contents

- Lead
- Guide wire
- Needles (epidural, modified Tuohy):
 - curved tip, 14 gauge, 8.9 cm
 - straight tip, 14 gauge, 8.9 cm
- Stylets:
 - bent, 0.30 mm diameter (green stylet handle and blue cap; inserted in the lead)
 - straight, 0.30 mm diameter (green stylet handle and green cap)
 - bent, 0.25 mm diameter (white stylet handle and blue cap)
 - straight, 0.25 mm diameter (white stylet handle and white cap)
- Product literature

Device specifications

The lead has electrodes on the distal end; the proximal (connector) end fits into an 8-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead.

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Table 1. Device specifications for leads: Models 977D160 and 977D260^a

Description	Models 977D160	Models 977D260
Connector	Octapolar, in-line	Octapolar, in-line
Conductor resistance ^{b,c}	1.6 Ω per cm	1.6 Ω per cm
Length	60 cm	60 cm
Diameter (lead body)	1.3 mm	1.3 mm
Distal (electrode) end		
Number of electrodes	8	8
Electrode shape	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm
Electrode spacing	1.5 mm	4.0 mm
Electrode distance	34.5 mm	52.0 mm
Distal tip distance	Not applicable	Not applicable
Radiopaque tip marker	No	No
Proximal (connector) end		
Lead contact length	1.5 mm	1.5 mm
'0' contact length	3.0 mm	3.0 mm
Lead contact spacing	2.8 mm	2.8 mm
Lead contact distance	22.5 mm	22.5 mm
Stylet handle length	29.0 mm	29.0 mm
Stylet wire diameter	0.25, 0.30 mm	0.25, 0.30 mm

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Table 2. Material of components in the Models 977D160 and 977D260 lead package

Component	Material	Material contacts human tissue
Lead		
Conductor wire	MP35N	No ^a
Conductor wire insulation	Fluoropolymer	No ^a
Electrodes	Platinum-iridium	Yes
Insulation	Polyurethane	Yes
Braiding	Tantalum	No ^a
Proximal contacts	MP35N	Yes
Lead tip	Epoxy	Yes
Adhesive	Epoxy	Yes
Guide wire	Stainless steel	Yes
Needles (curved, straight)	Stainless steel	Yes
Stylets (bent, straight)		
Handle	Polypropylene	Yes
Wire	Stainless steel	Yes

^a The material for this component has been tested for biocompatibility; however, the component is not intended to be in direct contact with human tissue.

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Cautions:

- Do not bend, kink, or stretch the lead or extension, which may damage the component.
- Do not use any instrument to handle the lead. The force may compress the lead, resulting in the inability to insert or withdraw the stylet.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

Cautions:

- Do not bend or kink the stylet, which may damage the component.
- When reinserting a stylet, use care because too much pressure on the stylet could damage the lead, resulting in intermittent or loss of stimulation.

Preparing for surgery

Before opening the lead package, verify the model number, use-by date, lead length, and connector type.

Placing a percutaneous lead

Warning: As with any spinal procedure, the risk of serious injury to the patient (eg, hemorrhage, hematoma, or paralysis) increases as location of the selected needle insertion site progresses up the vertebral column—from a lower risk at a lumbar location to a higher risk at a cervical location. Select a vertebral location that provides the widest and easiest access to the epidural space during needle insertion to reduce the risk of serious patient injury resulting from direct trauma to the spinal cord.

Caution: To reduce the risk of lead damage that may result in intermittent stimulation or loss of stimulation, requiring additional surgery to replace the lead:

- use only the modified Tuohy needle supplied in the kit.
 - use a shallow needle-insertion angle (45° or less) when inserting or withdrawing the needle into or out of the epidural space.
1. Make an incision at the needle-entry site to the depth of the subcutaneous fascia.
 2. Using a paramedial approach and under fluoroscopy, insert the needle (included in the kit) into the epidural space at the appropriate angle until you encounter resistance from the ligamentum flavum.
Note: Midline placement may cause passive damage to the lead over time because of ligament or spinous process movement.
 3. Confirm needle location under fluoroscopy.

4. After rotating the needle so that the beveled edge faces cephalad, remove the needle stylet.
5. Advance the needle and confirm entry into the epidural space (eg, using the loss-of-resistance technique with air or sterile [United States Pharmacopeia—USP] water).

Caution: Do not use contrast media or a saline flush. Contrast media may obscure the field of view and a saline flush may increase the difficulty of lead placement.

6. For a second lead, repeat steps 1 to 5 noting these recommendations:
 - Implant the second lead parallel to the first lead and approximately 1 to 3 mm lateral of the physiological midline.
 - Introduce the second lead one vertebral space below the first lead to help prevent nicking or cutting the first lead.
 - Stagger the lead tips or place them several vertebral spaces apart, depending on the position that produces the most effective paresthesia.

7. After inserting the guide wire through the needle, advance the guide wire no farther than 1 to 3 cm past the needle tip. Next, remove the guide wire from the needle.

Note: If the guide wire track deviates from the intended pathway, steering and manipulating the lead will be more difficult.

8. Using fluoroscopy, slowly insert the lead through the needle and advance the lead to the initial target placement site. A stylet may need to be reinserted.

Notes:

- The bent stylet is keyed so that the flat notch of the stylet handle faces the same direction as the distal bent tip.
- Ensure the lead electrodes are completely out of the distal end of the needle when changing a stylet. If the lead electrodes still partially remain in the needle, the force needed to insert the stylet may increase.
- The force to insert and withdraw the stylet increases as the needle angle increases.
- Before advancing the lead, ensure the stylet is fully inserted.
- If resistance is felt while steering (with the stylet inserted), rotate the lead and stylet or pull back slightly to ensure the lead is not stuck. Use short, firm movements to advance the lead.

9. After verifying the lead position under fluoroscopy (anterior-posterior and lateral views), compare that location with the location that has the highest probability of paresthesia coverage.

Notes:

- To increase lead stability, insert enough lead length to extend at least three vertebral bodies into the epidural space.

- Position the lead so that the center bipole electrodes (eg, electrodes 3 and 4 on an octapolar lead) will be active. If lead migration occurs, effective stimulation may be regained through programming adjustments instead of surgical repositioning.

Testing stimulation intraoperatively

1. Refer to the instructions associated with the external device(s) used in intraoperative test stimulation.
2. Identify optimal stimulation parameters, beginning at a pulse width of 210 to 240 μ s and a rate of 30 Hz.

Note: Ensure that the patient can provide immediate feedback.

3. Increase the amplitude beginning from 0.0 V while asking the patient closed-ended questions to identify the perception threshold (the amplitude at which the patient first perceives paresthesia), the discomfort threshold (the amplitude at which paresthesia is beyond the patient's tolerance), and the paresthesia coverage.

Note: If good paresthesia coverage is not attained, change electrode settings before repositioning the lead to confirm the direction of lead movement.

4. If two leads were placed, repeat steps 2 – 3 for the second lead; optimize paresthesia coverage using both leads.
5. In the patient's chart, document the lead position that provided appropriate stimulation coverage (ie, record the settings and patient responses and include a fluoroscopic image of the final lead position).

Removing the needle and closing

After the lead is in place, the physician can remove the needle, close, and use medical tape to secure the lead to the skin.

1. While holding the lead in place, disconnect the stylet from the lead. Partially expose, but do not completely withdraw the stylet.
2. While maintaining lead position, use minimal force to remove the needle and stylet.

Caution: Use minimal traction to remove the needle because quick or sudden removal may dislodge the lead.

3. If two leads were placed, repeat steps 1 – 2 for the second lead.
4. Connect the external device(s) used for trial screening.
Refer to the instructions associated with the external device(s) for connection information.
5. Verify test stimulation parameters to ensure that the lead has not moved. If the lead has moved, reposition it.
6. Close the lead exit site using standard suture technique.
7. Apply a dressing over the exit site of the lead.

Note: An antimicrobial dressing may be used at the lead exit site.

8. Place a sterile gauze pad between the skin and the lead. Place the lead over the sterile pad lateral to midline.
9. Use medical tape to secure the lead to the skin.

Preparing components for trial screening

The patient's response to the system and the efficacy of the system should be evaluated in both a clinical and a home setting. Different electrode configurations should be evaluated at various parameter settings during the trial evaluation period.

1. Place a gauze bandage on the skin where the lead and external device(s) used for trial screening will be placed on the patient.
2. Use the settings from "Testing stimulation intraoperatively" on page 13 to program the external neurostimulator for trial screening.
3. Tape the lead and the external device(s) together to secure the connection.
4. Tape the entire assembly to the patient's skin, allowing for strain relief.
5. Proceed with the trial evaluation.

Warning: Physicians should not prescribe MRI for patients undergoing trial stimulation or who have any neurostimulation system components that are not fully implanted. Explant all trial stimulation components if an MRI scan is required. MRI has not been tested on trial stimulation components and may cause heating of the lead electrodes, resulting in tissue damage or serious patient injury.

Removing components after trial screening

1. Verify that the external neurostimulator is turned off.
2. Refer to the instructions associated with the external device(s) used for trial screening for information on removing.
3. Remove all tape from the length of the lead body.
4. Open the original lead incision.
5. Use your thumb and index finger to grasp the lead as close as possible to the lead exit site.
6. Carefully retract the lead from the epidural space.
7. If the trial screening was successful, proceed with the chosen therapy.
8. If the trial screening was unsuccessful, clean the exit site(s) and close.

Medtronic

Alleviating Pain · Restoring Health · Extending Life

Manufacturer

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
www.medtronic.com
Tel. +1-763-505-5000
Fax +1-763-505-1000

Authorized Representative **in the European Community**

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. +31-45-566-8000
Fax +31-45-566-8668

Europe/Africa/Middle East **Headquarters**

Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
www.medtronic.eu
Tel. +41-21-802-7000
Fax +41-21-802-7900

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

On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Specify SureScan MRI 5-6-5 and 2x8 Lead Kit" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "1x8 low impedance extension kit" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
STERILE EO

Sterilization: ethylene-oxide gas


Caution: Consult accompanying documents

Manufacturing date

Use by

Storage temperature

SN

Serial number

Conformité Européenne (European Conformity).
This symbol means that the device fully complies
with European Directive AIMD 90/385/EEC.

! USA

For USA audiences only

Table of contents

Device description 5

Package contents 5

Device specifications 5

Instructions for use 6

Positioning the patient 6

Removing the remaining segment of the in-line percutaneous extension 7

Creating a subcutaneous pocket for the neurostimulator 7

Tunneling the extension 8

Creating a pocket for the lead-extension connection 9

Connecting the lead to the extension 9

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USA **Refer to the clinical summary booklet packaged with the neurostimulation system for information on the clinical study results of the neurostimulation system and individualization of treatment.**

37081 2007-10 English 3

Device description

The Medtronic Model 37081 1 x 8 Extension is part of a neurostimulation system for pain therapy.

Package contents

- Extension
- Connector boots, 1 x 8 low-profile (2)
- Setscrews (2)
- Tunneling tools: extension passer and obturator
- Wrench, torque
- Product literature
- Warranty card (USA only)
- Registration form (USA only)

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Figure 1. Model 37081 extension.

Table 1. Device specifications^a for the Model 37081 Extension.

Description	Model 37081
Resistance ^{b,c}	Maximum 16.0 Ω for all lengths
Length	10 – 110 cm
Contact spacing	2.8 mm
Distal end (lead)	
Connector	Octapolar, in-line
Lead entrance diameter	1.3 mm
Outer diameter	4.7 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Diameter	1.3 mm

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Table 2. Material of components in the Model 37081 package.

Component	Material	Material contacts human tissue
Extension		
Conductor wire	Silver core MP35N	No
Conductor wire insulation	Fluoropolymer	No
Insulation	Polyurethane, Silicone rubber	Yes
Setscrew connector block	Stainless steel	Yes
Spring contacts	MP35N, Platinum-iridium	Yes
Connector boots	Silicone rubber	Yes
Setscrew	Titanium	Yes
Tunneling tools		
Extension passer	Stainless steel and polypropylene	Yes
Obturator	Polypropylene	Yes
Wrench, torque		
Handle	Polymer	Yes
Shaft	Stainless steel	Yes

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Notes:

- Blunt dissection may be needed to expose the lead-percutaneous extension connection.
- If the percutaneous components were not removed, locate the wire inside the incision and cut the wire. Pull the outside percutaneous components and discard the cable and wires.

Caution: When severing the percutaneous extension, use minimal traction on the extension because strong force may dislodge the lead.

2. While maintaining lead position, carefully remove the lead-percutaneous extension connection from the incision.
3. Disconnect the lead from the percutaneous extension:
 - a. Remove the ligature, then slip the connector boot off the connection.

Caution: Do not use sharp instruments near the lead. Nicking or cutting the insulation may cause loss of stimulation or component failure requiring surgical replacement.

- b. Loosen the setscrew by turning the wrench counterclockwise.
 - c. Gently separate the extension setscrew connector from the lead.

Caution: If resistance is felt while removing the lead from the extension connector, first, loosen (but do not remove) the setscrew to ensure the lead contacts are not engaged. Next, inspect the lead connector for damage (eg, flattening of the lead contacts, stretching of the lead), which may result in intermittent or loss of stimulation.

- d. Remove and discard the connector boot; if necessary, carefully slit the boot.
4. Discard the remaining segment of the percutaneous extension.
5. For a second percutaneous extension, repeat steps 1 – 4.

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2. Assemble the tunneling tool and bend it as necessary. Next, tunnel subcutaneously from the neurostimulator pocket until the tip of the obturator is exposed at the lead site.

Cautions:

- Do not bend the extension passer at an angle greater than 90° at any one bend. An angle greater than 90° may damage the passer and prevent the extension from advancing, or damage the extension, resulting in intermittent or loss of stimulation.
- Proceed slowly when the tunneling tool approaches the lead site. If excess force is used, the patient could experience additional trauma when resistance to tunneling suddenly ceases.

Notes:

- Deep tunneling is not desirable.
 - Avoid the thoracic ribs.
 - If the tunneling tool does not extend to the lead site, make an intermediate incision.
3. Remove the obturator from the tunneling tool.
 4. Slide the proximal (neurostimulator) end of the extension through the extension passer until it exits at the handle or until fully encompassed by the extension passer. Repeat this step for two extensions.
 5. Remove the extension passer from the tunnel.
 6. Discard the tunneling tool.

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Caution: Do not place the lead-extension connection on the paraspinal muscles or the spinous processes. The connection may erode through the skin, resulting in infection and requiring removal of the neurostimulation system.

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Cautions:

- Do not use saline or other ionic fluids at connections, which could result in a short circuit.
- Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. Slip the connector boot onto the lead body with the open end facing out, exposing the lead contacts.

Notes:

- Irrigate the connector boot with a nonionic antibiotic solution. Dry the interior of the connector boot.
- A clear and white radiopaque boot are provided for use in dual-lead systems.

2. Wipe the lead body and the extension setscrew connector junction with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution. Dry all connections.
3. Insert the lead fully into the extension setscrew connector.

Notes:

- Each lead contact must be aligned under each extension connector contact.
- During insertion, some resistance is typical.

4. Use the torque wrench to tighten the setscrew to complete the electrical circuit.

Cautions:

- To avoid overtightening, do not use a hex wrench to tighten the extension setscrew. Overtightening the extension setscrew may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.
- Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.

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△ **Cautions:**

- Do not tie ligatures around the lead or extension body, which may damage the insulation.
- Do not overtighten ligatures on the anchor or connector boot. Ligatures that are too tight may damage the component.
- Ensure that the end of the connector boot is secure to prevent fluid in the lead-extension connection, which may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

7. For a second lead and extension, repeat steps 2 – 6.

Note: If two extensions were implanted but only one lead will be connected, cover the unused extension with a 1 x 8 low-profile closed connector boot from an accessory kit, and suture.

8. Using minimal force, pull the extension at the neurostimulator pocket and feed the lead-extension connection into the lead-extension connector pocket.

△ **Cautions:**

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 3). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 3. Coiling excess lead.

9. Insert the coiled lead body into the pocket, under the connection, leaving as much slack as possible in the lead body between the anchor and the lead-extension connection.
10. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.

Contacts:

Asia: Medtronic International Ltd.
Tel. 02891-4068
Fax 02591-0313

Asia: Medtronic Asia Ltd.
Tel. (02)-548-1148
Fax (02)-518-4786

Australia: Medtronic Australasia Pty. Ltd.
Tel. +61-2-9857-9000
Fax +61-2-9878-5100
Toll-free 1-800-668-670

Austria: Medtronic Österreich GmbH
Tel. (01)-240440
Fax (01)-24044-100

Belgium: Medtronic Belgium S.A.
Tel. 02-456-0900
Fax 02-460-2667

Canada: Medtronic of Canada Ltd.
Tel. (1905)-826-6020
Fax (1905)-826-6620

Czech Republic: Medtronic Czechia s.r.o.
Tel. 2-965-795-80
Fax 2-965-795-89

Denmark: Medtronic Danmark A/S
Tel. 45-32-48-18-00
Fax 45-32-48-18-01

Finland: Medtronic Finland Oy/Ltd
Tel. (09)-755-2500
Fax (09)-755-25018

France: Medtronic France S.A.S.
Tel. 01-5538-1700
Fax 01-5538-1800

Germany: Medtronic GmbH,
Tel. (02159)-81490
Fax (02159)-8149100

Greece: Medtronic Hellas S.A.
Tel. 210-67-79-099
Fax 210-67-79-399

Hungary: Medtronic Hungária Kft.
Tel. 1-889-06-00
Fax 1-889-06-99

Ireland: Medtronic Ireland Ltd.
Tel. (01)-890-6522
Fax (01)-890-7220

Italy: Medtronic Italia SpA
Tel. 02-241371
Fax 02-241381

Italy: Medtronic Italia SpA
Tel. 06-328141
Fax 06-3215812

Japan: Medtronic Japan
Tel. 3-6430-2011
Fax 3-6430-7140

Latin America: Medtronic, Inc.
Tel. (1305)-500-9328
Fax (1786)-709-4244

The Netherlands: Medtronic B.V.
Tel. (045)-566-8000
Fax (045)-566-8668

Norway: Medtronic Norge AS
Tel. 67-10-32-00
Fax 67-10-32-10

Poland: Medtronic Poland Sp. z o.o.,
Tel. (022)-465-69-00
Fax (022)-465-69-17

Portugal: Medtronic Portugal Lda.
Tel. 21-724-5100
Fax 21-724-5199

Russia: Medtronic Russia
Tel. (8495) 580-7377
Fax (8495) 580-7378

Spain: Medtronic Ibérica S.A.
Tel. 91-625-0400
Fax 91-650-7410

Sweden: Medtronic AB
Tel. 08-568-585-00
Fax 08-568-585-01

Switzerland: Medtronic (Schweiz) AG
Tel. 031-868-0100
Fax 031-868-0199

U.K.: Medtronic U.K. Ltd.
Tel. 01923-212213
Fax 01923-241004

USA: Medtronic, Inc.
Tel. (1763)-505-5000
Fax (1763)-505-1000
Toll-free (1800)-328-0810

Medtronic

Alleviating Pain · Restoring Health · Extending Life

Manufacturer

Medtronic Inc.
710 Medtronic Parkway,
Minneapolis MN 55432-5604,
USA.
Internet: www.medtronic.com
Tel. +1-763-505-5000
Fax +1-763-505-1000

Medtronic E.C. Authorized Representative/Distributed by

Medtronic B.V.
Earl Bakkenstraat 10,
6422 PJ Heerlen,
The Netherlands
Tel. +31-45-566-8000
Fax +31-45-566-8668

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl
Route du Molliat 31,
Case Postale
CH - 1131 Tolochenaz,
Switzerland
Internet: www.medtronic.com
Tel. +41-21-802-7000
Fax +41-21-802-7900

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* M A 0 4 3 2 7 A 0 0 6 *

© Medtronic, Inc. 2007
All Rights Reserved
MA04327A006/Rev A

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
Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA

Custodian:

Ashley Luloff
Regulatory Affairs Specialist
Medtronic Incorporated

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

On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Dual Quadripolar Extension Kit (8-2-4)" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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Tunneling the extension 10

Creating a pocket for the lead-extension connection 11

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USA Refer to the clinical summary booklet packaged with the neurostimulator for information on the clinical study results of the neurostimulation system and individualization of treatment.

37082 2008-06 English 5

Description

The Medtronic Model 37082 Extension is part of a neurostimulation system for pain therapy

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
Figure 1. Model 37082 extension.

Note: The marker band is for lead identification. The lead connected to the extension leg with the marker band is for either circuits 0-3 or circuits 8-11.

Table 1. Device specifications^a for the Model 37082 Extension

Description	Model 37082
Resistance ^{b,c}	Maximum 24.0 Ω for all lengths
Length	10 – 110 cm
Distal end (lead)	
Connectors	Quadripolar, in-line (2)
Contact spacing	4.3 mm
Lead entrance diameter	1.5 mm
Connector outer diameter	3.8 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Contact spacing	2.8 mm
Diameter	1.3 mm

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^c Resistance is proportional to length; long lengths have higher resistance that may limit the amplitude.

Table 2. Material of components in the Model 37082 package

Component	Material	Material contacts human tissue
Extension		
Conductor wire	Silver core MP35N	No
Distal end (lead)		
Insulation	Silicone rubber	Yes
Marker band	Platinum-iridium	No
Setscrew connector blocks	Stainless steel	Yes
Proximal end (neurostimulator)		
Conductor wire insulation	Fluoropolymer	No
Contacts	MP35N	Yes
Insulation	Polyurethane	Yes
Connector boots	Silicone rubber	Yes
Neurostimulator plug	Polyurethane	Yes
Setscrews	Titanium	Yes
Tunneling tools		
Carrier	Polypropylene	Yes
Extension passer	Stainless steel and polypropylene	Yes
Obturator	Polypropylene	Yes
Wrench, torque		
Handle	Polymer	Yes
Shaft	Stainless steel	Yes

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 **Cautions:**

- When routing the extension, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When multiple leads are implanted, route the lead-extensions so the area between them is minimized (Figure 2). If the lead-extensions are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).


Figure 2. Routing for multiple lead-extensions.

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Notes:

- Deep tunneling is not desirable.
 - Avoid the thoracic ribs.
 - If the tunneling tool does not extend to the neurostimulator pocket, make an intermediate incision.
3. Remove the obturator from the tunneling tool.
 4. Insert the carrier into the extension passer.
 5. Insert the extension setscrew connector junctions into the carrier grooves. Insert the first extension into the carrier groove, then gently pull the proximal end of the carrier further through the extension passer handle and snap in place. Gently place the second extension in the other carrier groove.

△ Caution: Use care when inserting the extension body into the carrier groove. Rough handling may damage the extension insulation.

6. Using minimal force, pull the carrier-extension assembly through the tunnel to where the lead is anchored.
Note: Ensure that the bifurcation portion of the extension remains in the neurostimulator pocket when the extension is tunneled.
7. Remove the extension setscrew connector junctions from the carrier grooves.
8. Discard the tunneling tool.

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△ Cautions:

- To avoid overtightening, use only the torque wrench included in the extension package to tighten the extension setscrews. The torque wrench included in the extension package is marked with an icon to indicate it is not for use with the neurostimulator. Overtightening the extension setscrews may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.
- Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.

5. Using minimal force, pull on the connection to ensure that it is secure.
6. Slide the connector boot over the lead-extension connection and tie a ligature (using nonabsorbable silk) between the pairs of rings at either end. Suture the ends of the connector boot.

Cautions:

- Do not tie ligatures around the lead or extension body, which may damage the insulation.
- Do not overtighten ligatures on the connector boot. Ligatures that are too tight may damage the connector boot, the lead, or the extension.
- Ensure that the end of the connector boot is secure to prevent fluid in the lead-extension connection, which may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

7. For a second lead, repeat steps 2 – 6.

Note: If only one lead will be connected to the extension, cover the unused connector with a 1 x 4 low-profile closed connector boot from an accessory kit, and suture.

Note: When using one or two leads, use the marker band on the leg of the extension as a reference for lead identification.

8. Using minimal force, pull the extension at the neurostimulator pocket and feed the lead-extension connections into the lead-extension connector pocket.

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
Figure 4. Coiling excess lead.

9. Insert the coiled lead body into the pocket, under the connections, leaving as much slack as possible in the lead body between the anchor and the lead-extension connections.
10. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.

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Medtronic
Further Together

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA

Ashley Luloff
Custodian:
Ashley Luloff
Regulatory Affairs Specialist
Medtronic Incorporated

July 2, 2019


Instruction for use

For

Single Quadripolar Extension Kit (8-4)

ATTESTATION

STATE OF MINNESOTA
HENNEPIN COUNTY

On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Single Quadripolar Extension Kit (8-4)" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Pamela M. Heitz
Official Notary Signature

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Rx only

2005

Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

Extension length

Do not reuse

Sterilization: ethylene-oxide gas

Caution, consult accompanying documents


Manufacturing date

Manufacturer

Use by

Storage temperature

Serial number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC.

Authorized representative in the European community

For USA audiences only

Table of contents

Description 7

Package contents 7

Device specifications 7

Instructions for use 9

Positioning the patient 9

Removing the remaining segment of the in-line percutaneous extension 9

Creating a subcutaneous pocket for the neurostimulator 10

Tunneling the extension 10

Creating a pocket for the lead-extension connection 12

Connecting the lead to the extension 12

Refer to the indications sheet for indications and related information.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.

Refer to System Eligibility, Battery Longevity, Specifications reference manual packaged with the software application card for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.

 Refer to the clinical summary booklet packaged with the neurostimulator for information on the clinical study results of the neurostimulation system and individualization of treatment.

Description

The Medtronic Model 37083 Extension is part of a neurostimulation system for pain therapy.

Package contents

- Extension
- Connector boots, 1 x 4 low-profile (2)
- Octapolar in-line neurostimulator plug
- Setscrews (2)
- Tunneling tools: extension passer, obturator, carrier
- Wrench, torque
- Product literature
-  Warranty card
-  Registration form

Device specifications

The extension has four connectors on the distal end and eight connectors on the proximal end. The in-line distal end connects to a Medtronic lead. The proximal end connects to a Medtronic neurostimulator.


Figure 1. Model 37083 extension.

Table 1. Device specifications^a for the Model 37083 Extension

Description	Model 37083
Resistance ^{b,c}	Maximum 24.0 Ω for all lengths
Length	10 – 110 cm
Distal end (lead)	
Connector	Quadripolar, in-line
Contact spacing	4.3 mm
Lead entrance diameter	1.5 mm
Outer diameter	3.8 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Contact spacing	2.8 mm
Diameter	1.3 mm

^a All measurements are approximate.

^b Electrical resistance of this device only.

^c Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

Table 2. Material of components in the Model 37083 package

Component	Material	Material contacts human tissue
Extension		
Conductor wire	Silver core MP35N	No
Distal end (lead)		
Insulation	Silicone rubber	Yes
Setscrew connector block	Stainless steel	Yes
Proximal end (neurostimulator)		
Conductor wire insulation	Fluoropolymer	No
Contacts	MP35N	Yes
Insulation	Polyurethane	Yes
Connector boots	Silicone rubber	Yes
Neurostimulator plug	Polyurethane	Yes
Setscrews	Titanium	Yes
Tunneling tools		
Carrier	Polypropylene	Yes
Extension passer	Stainless steel and polypropylene	Yes
Obturator	Polypropylene	Yes
Wrench, torque		
Handle	Polymer	Yes
Shaft	Stainless steel	Yes

Instructions for use

Implanting physicians should be experienced in epidural-access procedures and should be thoroughly familiar with all product labeling.

When the patient completes a test stimulation period to determine the effectiveness of stimulation in managing pain, the percutaneous components used to connect the snap-lid connector cable to the lead may be removed at the follow-up visit.

Cautions:

- Do not bend, kink, or stretch the lead or extension body, which may damage the component.
- Do not use any instrument to handle the extension. The force may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.

Positioning the patient

1. Select a location for the neurostimulator and position the patient accordingly.

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Note: Also take into account:

- possible future cardiac needs (eg, pacemaker, defibrillator). Implanting a neurostimulator on the patient's right side allows for possible future placement of cardiac devices on the patient's left side.
 - the cosmetic needs of the patient.
2. Before opening the extension package, verify the model number, use-by date, extension length, and connector type.
 3. Identify the neurostimulator pocket site and mark the intended extension route from the lead-incision site to the neurostimulator pocket.

Removing the remaining segment of the in-line percutaneous extension

1. After finding and exposing the percutaneous extension connector, make an incision large enough to grasp the lead body and maintain lead position.

Notes:

- Blunt dissection may be needed to expose the lead-extension connection.
- If the percutaneous components were not removed, locate the wire inside the incision and cut the wire. Pull the outside percutaneous components and discard the cable and wires.

 **Caution:** When severing the percutaneous extension, use minimal traction on the extension because strong force may dislodge the lead.

2. While maintaining lead position, carefully remove the lead-extension connection from the incision.
3. Disconnect the lead from the percutaneous extension:
 - a. Remove the ligature, then slip the connector boot off the connection.

 **Caution:** Do not use sharp instruments near the lead. Nicking or cutting the insulation may cause loss of stimulation or component failure requiring surgical replacement.

- b. Loosen each setscrew by turning the wrench counterclockwise.
 - c. Gently separate the extension setscrew connector from the lead.

 **Caution:** If resistance is felt while removing the lead from the extension connector, first, loosen (but do not remove) the setscrews to ensure the lead contacts are not engaged. Next, inspect the lead connector for damage (eg, flattening of the lead contacts, stretching of the lead), which may result in intermittent or loss of stimulation.

- d. Remove and discard the connector boot; if necessary, carefully slit the boot along the extension setscrews.
 4. Discard the remaining segment of the percutaneous extension.
 5. For a second percutaneous extension, repeat steps 1–4.

Creating a subcutaneous pocket for the neurostimulator

1. Make an incision equal to the length of the neurostimulator.
2. Use blunt dissection to create a subcutaneous pocket.

Note: Refer to the neurostimulator implant manual for the proper subcutaneous pocket depth.

Tunneling the extension

For single extension procedures, tunneling can be performed without or with a carrier. For dual extension procedures, the carrier must be used.

1. Identify the tunneling route between the lead incision and the neurostimulator pocket.

 **Cautions:**

- When routing the extension, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When two leads are implanted, route the lead-extensions so the area between them is minimized (Figure 2). If the lead-extensions are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).

Proper

Improper


Figure 2. Routing for two lead-extensions.

2. Assemble the tunneling tool and bend it as necessary. Next, tunnel subcutaneously from the lead site to the neurostimulator site until the tip of the obturator is exposed at the neurostimulator pocket.

Cautions:

- Do not bend the extension passer at an angle greater than 90° at any one bend. An angle greater than 90° may damage the passer and prevent the extension from advancing, or damage the extension, resulting in intermittent or loss of stimulation.
- Proceed slowly when the tunneling tool approaches the neurostimulator pocket. If excess force is used, the patient could experience additional trauma when resistance to tunneling suddenly ceases.

Notes:

- Deep tunneling is not desirable.
- Avoid the thoracic ribs.
- If the tunneling tool does not extend to the neurostimulator pocket, make an intermediate incision.

3. Remove the obturator from the tunneling tool.
4. If using the carrier, insert it into the extension passer.

Note: For a single extension procedure, snap the proximal end of the carrier into the extension passer handle to secure it in place.

5. Insert the extension setscrew connector junction(s) into the carrier groove(s) or the extension passer.

If the carrier is not used – slide the distal end of the extension from the neurostimulator pocket towards the handle of the extension passer while keeping enough of the proximal end (neurostimulator connector end) of the extension available to grip during the removal of the passer.

For dual extension procedures – Insert the first extension into the carrier groove, then gently pull the proximal end of the carrier further through the extension passer handle and snap in place. Gently place the second extension in the other carrier groove.

Caution: Use care when inserting the extension body into the carrier groove. Rough handling may damage the extension insulation.

6. Using minimal force, pull the carrier-extension assembly through the tunnel to where the lead is anchored. If the carrier is not used, simply remove the extension passer from the tunnel.

Note: Ensure that the transition portion of the extension remains in the neurostimulator pocket when the extension is tunneled.

7. As required, remove the extension setscrew connector junction(s) from the carrier groove(s).
8. Discard the tunneling tool.

Creating a pocket for the lead-extension connection

Use blunt dissection to form a subcutaneous lead-extension connector pocket. The pocket should be:

- lateral to the lead incision.
- deep enough to minimize the potential for skin erosion or irritation.
- large enough to accommodate the lead-extension connection(s).

△ **Caution:** Do not place the lead-extension connection on the paraspinal muscles or the spinous processes. The connection may erode through the skin, resulting in infection and requiring removal of the neurostimulation system.

Connecting the lead to the extension

△ **Cautions:**

- Do not use saline or other ionic fluids at connections, which could result in a short circuit.
 - Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.
1. Slip the connector boot onto the lead body with the open end facing out, exposing the lead contacts.

Notes:

- Irrigate the connector boot with a nonionic antibiotic solution. Dry the interior of the connector boot.
 - A clear and white radiopaque boot are provided for use in dual-lead systems.
2. Wipe the lead body and the extension setscrew connector junction with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution. Dry all connections.
 3. Insert the lead fully into the extension setscrew connector.
- Note:** Each lead contact must be aligned under each setscrew.
4. Use the torque wrench to tighten each setscrew to complete the electrical circuit. Hold the extension firmly on the sides of each connector block while tightening the setscrew (Figure 3).


Figure 3. Tighten the setscrews holding the extension firmly on the sides of each connector block.

△ **Cautions:**

- To avoid overtightening, use only the torque wrench included in the extension package to tighten the extension setscrews.

The torque wrench included in the extension package is marked with an icon to indicate it is not for use with the neurostimulator. Overtightening the extension setscrews may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.

- Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.
5. Using minimal force, pull on the connection to ensure that it is secure.
 6. Slide the connector boot over the lead-extension connection and tie a ligature (using nonabsorbable silk) between the pairs of rings at either end. Suture the ends of the connector boot.

Cautions:

- Do not tie ligatures around the lead or extension body, which may damage the insulation.
 - Do not overtighten ligatures on the connector boot. Ligatures that are too tight may damage the connector boot, the lead, or the extension.
 - Ensure that the end of the connector boot is secure to prevent fluid in the lead-extension connection, which may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.
7. For a second lead and extension, repeat steps 2 – 6.
- Note:** If two extensions were implanted but only one lead will be connected, cover the unused extension with a 1 x 4 low-profile closed connector boot from an accessory kit, and suture.
8. Using minimal force, pull the extension at the neurostimulator pocket and feed the lead-extension connection into the lead-extension connector pocket.

Cautions:

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 4). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 4. Coiling excess lead.

9. Insert the coiled lead body into the pocket, under the connection, leaving as much slack as possible in the lead body between the anchor and the lead-extension connection.

10. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.

Contacts:**Asia:**

Medtronic International Ltd.
Tel. 02919-1362
Fax 02907-3998

Medtronic Asia Ltd.

Tel. (02)-548-1148
Fax (02)-518-4786

Australia:

Medtronic Australasia Pty. Ltd.
97 Waterloo Road
North Ryde, NSW 2113
Australia
Tel. +61-2-9857-9000
Fax +61-2-9878-5100
Toll free 1-800-668-6700

Austria:

Medtronic Österreich GmbH
Tel. 01-240440
Fax 01-24044-100

Belgium:

Medtronic Belgium S.A.
Tel. 02-456-0900
Fax 02-460-2667

Canada:

Medtronic of Canada Ltd.
Tel. (1-905)-460-3800
Fax (1905)-826-6620

Czech Republic:

Medtronic Czechia s.r.o.
Tel. 2-965-795-80
Fax 2-965-795-89

Denmark:

Medtronic Danmark A/S
Tel. 45-32-48-18-00
Fax 45-32-48-18-01

Finland:

Medtronic Finland Oy/LTD
Tel. (09)-755-2500
Fax (09)-755-25018

France:

Medtronic France S.A.S.
Tel. 01-5538-1700
Fax 01-5538-1800

Germany:

Medtronic GmbH
Tel. (02159)-81490
Fax (02159)-8149100

Greece:

Medtronic Hellas S.A.
Tel. 210-67-79-099
Fax 210-67-79-399

Hungary:

Medtronic Hungária Kft.
Tel. 1-889-06-00
Fax 1-889-06-99

Ireland:

Medtronic Ireland Ltd.
Tel. (01)-890-6522
Fax (01)-890-7220

Italy:

Medtronic Italia SpA
Tel. 02-241371
Fax 02-241381

Tel. 06-328141

Fax 06-3215812

Japan:

Medtronic Japan
Tel. 03-6430-2016
Fax 03-6430-7110

Latin America:

Medtronic, Inc.
Tel. (1305)-500-9328
Fax (1786)-709-4244

Norway:

Medtronic Norge AS
Tel. 067-10-32-00
Fax 067-10-32-10

Poland:

Medtronic Poland Sp. z o.o.
Tel. (022)-465-69-00
Fax (022)-465-69-17

Portugal:

Medtronic Portugal, Lda.
Tel. 21-724-5100
Fax 21-724-5199

Russia:

Medtronic Russia
Tel. (8495) 580-7377
Fax (8495) 580-7378

Slovakia

Medtronic Slovakia, o.z.
Tel. 0268 206 911
Fax 0268 206 999

Spain:

Medtronic Ibérica, S.A.
Tel. 91-625-0400
Fax 91-650-7410

Sweden:

Medtronic AB
Tel. 08-568-585-00
Fax 08-568-585-01

Switzerland:

Medtronic (Schweiz) AG
Tel. 031-868-0100
Fax 031-868-0199

The Netherlands:

Medtronic B.V.
Tel. (045)-566-8000
Fax (045)-566-8668

U.K.:

Medtronic U.K. Ltd.
Tel. 01923-212213
Fax 01923-241004

USA:

Medtronic, Inc.
Tel. (1763)-505-5000
Fax (1763)-505-1000
Toll-free: (1-800)-328-0810

Medtronic

Alleviating Pain · Restoring Health · Extending Life

Manufacturer

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
www.medtronic.com
Tel. 1-763-505-5000
Fax 1-763-505-1000

Authorized Representative
in the European Community

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Europe/Africa/Middle East

Headquarters

Medtronic International Trading Sàrl
Route du Molliat 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
www.medtronic.eu
Tel. 41-21-802-7000
Fax 41-21-802-7900

Asia-Pacific

Medtronic International Ltd.
Suite 1106-11, 11/F, Tower 1, The Gateway,
25 Canton Road, Tsimshatsui,
Kowloon,
Hong Kong
Tel. +852-2919-1362
Fax +852-2907-3998

Contacts for specific countries are listed inside this cover.

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
Medtronic
Further. Together

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA

Custodian
Ashley Luloff
Regulatory Affairs Specialist
Medtronic Incorporated

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

On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "3550 3555 Accessory Kits" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

3550
3555

1998 (3550-01, -02
-03, -04, -05, -06,
-07, -08, -09)
2000 (3555)
2002 (3550-17, -22
-23, -24, -25)

**Medtronic E.C. Authorized
Representative/Distributed By**
Medtronic B.V.
Earl Bekkenstraat 10,
6422 PJ Heerlen,
The Netherlands.
Tel. +31-45-566-8000.
Fax +31-45-566-8668.

EC	REP
----	-----

MA00243A004 Rev. A
© Medtronic, Inc. 2011
All Rights Reserved
February 2011

MA00243A004 Rev. A

[illegible]

Kit doesn't contain this sheet; shown for reference only!

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
On this 8th day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Multi-Lead Trialing Cable Model 355531" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

July 8, 2019

Multi-lead Trialing Cable

355531

Trialing and intraoperative screening cable

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
Medtronic
Further Together

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA

Custodian:
Ashley Luloff
Regulatory Affairs Specialist
Medtronic Incorporated

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
On this 8th day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Passing Elevator Accessory Kit" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "3550-26 and 3550-27 Accessory Kits" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

Accessory Kits 3550-26
Kits d'accessoires 3550-27
Zubehörssets
Equipos de accesorios
Accessory kits
Kit degli accessori
Tilbehørssett
Tilbehørsset
Σετ παρελκομένων
Tilbehørsset
Conjuntos de accesorios
Zestawy akcesoriów
Sady príslušenství
Kiegészítő kiegészítők
Комплекты аксессуаров
Supravy príslušenstva
Aksesuar Setleri
Kituri de accesorii

Referral sheet * Feuille de référence
Referenzblatt * Hoja de referencia * Referenzblatt
Foglio di riferimento * Referenzark
Referenzblatt * Evropského evropského
Hemingshøjde * Fiches de référence * Symbole
adresy * Referenční arch * Információs lap
Лист напаванья * Referenční nárok * Базуру
seyfisi * Fyß de referenji

USA Rx only

Manufacturer
Medtronic Inc
710 Medtronic Parkway
Minneapolis MN 55432-5604
USA
www.medtronic.com
Tel: +1-763-505-6000
Fax: +1-763-505-1000

Authorized Representative
in the European Community
Medtronic B.V.
Een Bakkenstraat 10,
6422 PJ Heerlen,
The Netherlands
Tel: +31-45-566-8000
Fax: +31-45-566-8668

Europe/Africa/Middle East
Headquarters
Medtronic International Trading Sàrl
Route du Molin 31
Case Postale 84
CH - 1131 Tolochenaz,
Switzerland
www.medtronic.eu
Tel: +41-21-602-7000
Fax: +41-21-602-7900

Asia-Pacific
Medtronic International Ltd
Suite 1106-11, 11/F, Tower 1, The Gateway,
25 Canton Road, Tsimshatsui,
Kowloon
Hong Kong
Tel: +852-2919-1362
Fax: +852-2907-3998

© Medtronic, Inc. 2012
All Rights Reserved
MA13108A007

English

Models 3550-26 and 3550-27 1x8 Percutaneous Lead Revision Accessory Kits

The Models 3550-26 and 3550-27 1x8 Percutaneous Lead Revision Accessory Kits are intended for use when revising or repositioning an implanted Medtronic Neurological 1x8 percutaneous lead. With the exception of the lead and percutaneous extension, this kit contains all the sterile accessories supplied in lead kits for the following models:

- 3776 / 3876 1x8 SC
- 3777 / 3877 1x8
- 3778 / 3878 1x8 Compact

For information on the use of the accessories, refer to the lead implant manual included in this package.

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Instruction for use


For

Percutaneous Extension and Tunneling Tool Accessory Kit for Deep Brain Stimulation

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On this 8th day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Percutaneous Extension and Tunneling Tool Accessory Kit for Deep Brain Stimulation" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

July 8, 2019

Percutaneous Extension and Tunneling Tools Accessory Kit
For Deep Brain Stimulation

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Instruction for use


For

Curved Tip Epidural Needle Accessory Kits

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On this 8th day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Curved Tip Epidural Needle Accessory Kits" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

July 8, 2019

Medtronic

Curved-tip epidural needle 3550-14
accessory kits 3550-43

Kits d'accessoires d'aiguilles épidurales avec extrémité courbée
Epiduralknüle-Zubehörsatz (gebogene Spitzen)
Kit de accesorios de aguja epidural de punta curva
Toesbehorensets met gebogen epiduraalnaald
Kit degli accessori per ago epidurale a punta ricurva
Tilbehørssett for epiduralkanyler med bøyd spiss
Tilbehörsset för epiduralnålar med böjd spets
Κη παρεκδομήτων βελόνων επικληρίδιου με κυρτό άκρο
Tilbehørssett med epiduralkanyler med bøjet spids
Conjuntos de acessórios de agulhas epidurais com ponta curva
Hajlított végű epidurális tű kiegészítő készletei
Zestawy akcesoriów igieł nadoponowych o zakrzywionych końcówkach
Sady príslušenství epidurálních jehel s ohnutým hrotem
Комплекты принадлежностей эпидуральной иглы с изогнутым кончиком
Súpravy príslušenstva k epidurálnym ihlám zo zakriveným hrotom
Eğri uçlu epidural iğne aksesuar setleri
Kituri de accesorii cu ace epidurale cu vârf curbă

Referral sheet • Feuille de référence • Referenzblatt • Hoja de referencia • Referenzblatt • Foglietto di riferimento • Referansark • Referensblad • Ενημερωτικό έντυπο • Henvisningsblad • Folheto de referência • Információs lap • Informacje o dokumentacji • Referenční listák • Информационный бюллетень • Referenčný hárok • Başvuru sayfası • Fiş de recomandări

2009

Medtronic

Advancing Care • Restoring Health • Extending Life

Manufacturer
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Internet: www.medtronic.com
Tel. 1-763-505-5000
Fax 1-763-505-1000

Medtronic E.C. Authorized Representative/Distributed by
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Europe/Africa/Middle East Headquarters
Medtronic International Trading Sàrl
Route du Molliat 31
Case Postale
CH-1131 Tolochenaz
Switzerland
Internet: www.medtronic.co.uk
Tel. 41-21-802-7000
Fax 41-21-802-7900

Asia-Pacific
Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2891-4068
Fax 852-2591-0313

* 9 3 0 0 8 1 A 0 0 3 *

© Medtronic, Inc. 2009
All Rights Reserved
M930081A003

Explanation of symbols on package labeling. Refer to the appropriate product for symbols that apply.

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Description

The 14-gauge epidural needles contained in the Medtronic™ Neurostimulation accessory kit models 3550-14 and 3550-43 are designed for use with Medtronic percutaneous leads for spinal cord stimulation.

For information and instructions on the use of these accessories and tools including indications, contraindications, warnings, precautions, and adverse events, refer to the appropriate lead implant manual.

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
Custodian:
Ashley Luloff

Regulatory Affairs Specialist
Medtronic Incorporated

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
On this 8th day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Accessory Kits, Models 9779060,97790975,9779090" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

July 8, 2019

Accessory Kits
9779060
9779075
9779090

Kits d'accessoires
Zubehörsets
Kits de accesorios
Toebehorenssets
Kit degli accessori
Tilbehørssett
Tillbehörssset
Κιτ παρελκομένων
Tilbehørssætene
Conjuntos de acessórios
Zestawy akcesoriów
Sady příslušenství
Tartozékkészletek
Наборы принадлежностей
Súpravy s príslušenstvom
Aksesuar Setleri
Kifurli de accesorii

Referral sheet • Feuille de référence •
Referenzblatt • Hoja de referencia • Referentiblad •
Foglioletto di riferimento • Referanseark •
Referensblad • Φύλλο παραπομπής •
Henvisningsblad • Folheto de referência •
Informacje o dokumentacji • Referenční listák •
Információs lap • Информационный бюллетень •
Referentný hárok • Yönlendirme sayfası • Fişâ de
referință

USA Rx only

2013

Manufacturer
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5804
USA
www.medtronic.com
Tel. +1-763-505-5000
Fax +1-763-505-1000

Authorized Representative 
In the European Community
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. +31-45-566-8000
Fax +31-45-566-8668

**Europe/Africa/Middle East
Headquarters**
Medtronic International Trading Sàrl
Route du Molliat 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
www.medtronic.eu
Tel. +41-21-802-7000
Fax +41-21-802-7900

Asia-Pacific
Medtronic International Ltd.
Suite 1106-11, 11/F, Tower 1, The Gateway
25 Canton Road, Tsimshatsui
Kowloon
Hong Kong
Tel. +852-2919-1362
Fax +852-2907-3998

© Medtronic, Inc. 2012
All Rights Reserved
M939861A002

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Description

Models 9779060, 9779075, and 9779090
Revision Accessory Kits for the Vectris™
SureScan™ MRI 1x8 Percutaneous Leads

The Models 9779060, 9779075, and 9779090 Revision Accessory Kits are intended for use when revising or repositioning an implanted Medtronic Neuromodulation Vectris™ SureScan™ MRI 1x8 percutaneous lead. With the exception of the lead, each kit contains all the sterile accessories supplied in the Vectris SureScan MRI lead kits according to lead length.

The Model 9779060 Revision Accessory Kit is for use with the following Vectris SureScan MRI lead models (60 cm in length):

- 977A160 1x8 Subcompact
- 977A260 1x8 Compact

The Model 9779075 Revision Accessory Kit is for use with the following Vectris SureScan MRI lead models (75 cm in length):

- 977A175 1x8 Subcompact
- 977A275 1x8 Compact

The Model 9779090 Revision Accessory Kit is for use with the following Vectris SureScan MRI lead models (90 cm in length):

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On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Injex Bumpy Anchor" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
On this 2nd day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Injex Bi-Wing Anchor" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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

On this 8th day of July, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following "Accessory Kits" Instruction for use, issued by Medtronic Inc., presented to me by the document's custodian, Ashley Luloff, and, to the best of my knowledge, that the photocopies document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

July 8, 2019

Medtronic

Accessory Kits
Kits d'accessoires
Zubehörsets
Equipo de accesorios
Accessory kits
Kit degli accessori
Tilbehørssett
Tilbehørssett
Σετ παρελκόμενων
Tilbehørsæt
Zestawy akcesoriów
Conjuntos de acessórios
Комплекты принадлежностей
Aksesuar Kitleri
Soupravy příslušenství
Tartozékkészletek
Súpravy príslušenstva
Kituri de accesorii

Referral Sheet • Feuille de référence
Referenzblatt • Hoja de referencia
Referentiesheet • Foglio di riferimento
Referansark • Referensblad • Παραπεμπτικό έντυπο
Henvisningsblad • Ulotka informacyjna
Folheto de referência
Информационный бюллетень
Başvuru Sayfası • Referenční list • Adatlap
Referenčný hárok • Referral Sheet

Medtronic
Advancing Life - Saving Health - Transforming Life

Manufacturer
Medtronic, Inc.
710 Medtronic Parkway,
Minneapolis MN 55432-5604,
USA
Internet: www.medtronic.com
Tel. +1-763-505-5000
Fax +1-763-505-1000

Medtronic E.C.
Authorized Representative/
Distributed By
Medtronic B.V.
Earl Bakkenstraat 10,
6422 PJ Heerlen,
The Netherlands
Tel. +31-45-566-8000
Fax +31-45-566-8668

Europe/Africa/Middle East Headquarters
Medtronic International Trading Sarl
Route du Molliat 31,
Case Postale
CH - 1191 Tolochenaz,
Switzerland
Internet: www.medtronic.co.uk
Tel. +41-21-802-7000
Fax +41-21-802-7900

Asia-Pacific
Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza,
The Lee Gardens, 33 Hysan Avenue,
Causeway Bay,
Hong Kong
Tel. +852-2591-4068
Fax +852-2591-0313

© Medtronic, Inc. 2008
All Rights Reserved
MA08933A005/Rev A
* MA 08933A 005 *

English

Neurostimulation Accessory Kits

The components contained in the Medtronic Neurostimulation Accessory Kits are spare parts. For information and instructions on the use of these accessories and tools including indications, contraindications, warnings, precautions, and adverse events, refer to the appropriate product manual.

Important: When replacing a component with a spare part from an accessory kit, make sure the component in the accessory kit matches the original component.

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