Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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Important information for users of KARL STORZ devices

1 Important information for users of KARL STORZ devices

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This instruction manual is intended to aid you in the proper installation, connection, and operation of this device. All required details and all actions on your part are clearly presented and explained. Please read this manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

NOTE: This IFU is applicable for IMAGE1 S CONNECT® TC 200 and for IMAGE1 S CONNECT® II TC 201. Only when a special feature of IMAGE1 S CONNECT® TC 200 or IMAGE1 S CONNECT® II TC 201 is described the model number TC 200 or TC 201 is mentioned in the IFU to distinguish between CONNECT and CONNECT II. Otherwise the term «IMAGE1 S CONNECT®» is used.

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2 Images of the equipment

2.1 Possible combinations for IMAGE1 S™

The IMAGE1 S™ comprises IMAGE1 S CONNECT® (TC 200) or IMAGE1 S CONNECT® II (TC 201) and a minimum of one and maximum of three additional LINK modules. The IMAGE1 S™ H3-LINK (TC 300), IMAGE1 S™ X-LINK (TC 301), IMAGE1 S D3-LINK™ module (TC 302) and IMAGE1 S™ 4U-LINK module (TC 304) are available as LINK modules.

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Controls, displays, connectors, and their uses

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3 Controls, displays, connectors, and their uses

3.1 IMAGE1 S CONNECT® TC 200 in combination with the H3-LINK, X-LINK and D3-LINK modules (see page IV)

- ① Power switch
- ② White balance
- ③ USB port
- ④ Potential equalization connector
- ⑤ LINK inputs
- ⑥ SCB input/output sockets
(KARL STORZ Communication Bus)
- ⑦ LAN Interface
- ⑧ USB interface
- ⑨ HD DVI-D video output sockets
- ⑩ Potential equalization connector
- ⑪ 3G-SDI video output
- ⑫ Power cord socket
- ⑬ LINK output
- ⑭ H3-LINK camera connector
- ⑮ X-LINK camera connector
- ⑯ D3-LINK camera connector

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**3.2 IMAGE1 S CONNECT® II
TC 201 in combination with the
4U-LINK, X-LINK and D3-LINK
modules (see page V)**

- ① Power switch
- ② White balance
- ③ USB port
- ④ Potential equalization connector
- ⑤ LINK inputs
- ⑥ SCB input/output socket
(KARL STORZ Communication Bus)
- ⑦ LAN interface
- ⑧ **USB** interface
- ⑨ HD DVI-D video output socket
- ⑩ Potential equalization connector
- ⑪ 12G/3G-SDI video output
- ⑫ Power cord socket
- ⑬ LINK output
- ⑭ 4U-LINK camera connector
- ⑮ X-LINK camera connector
- ⑯ D3-LINK camera connector
- ⑰ DisplayPort video output sockets
- ⑱ EXT SYNC connector (for future use)

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Non-ionizing electromagnetic radiation

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This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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1	Important information for users of KARL STORZ devices	III
2	Images of the equipment	IV
2.1	Possible combinations for IMAGE1 S™	IV
3	Controls, displays, connectors, and their uses	VI
3.1	IMAGE1 S CONNECT® TC 200 in combination with the H3-LINK, X-LINK and D3-LINK modules (see page IV)	VI
3.2	IMAGE1 S CONNECT® II TC 201 in combination with the 4U-LINK, X-LINK and D3-LINK modules (see page V)	VII
4	Symbols employed	VIII
4.1	Symbols for operation	VIII
4.2	Label and packaging symbols	IX
5	General information	4
5.1	Description of the devices	4
5.2	Property rights	4
6	Safety instructions	5
6.1	Explanation of warnings and cautions	5
6.2	General safety information	6
6.3	Indications for use	10
6.3.1	Intended use	10
6.4	Contraindications	10
6.5	Patient profile	11
6.6	User profile (specialist physician)	11
6.7	User profile (assistant)	12
6.8	Safety precautions at the site of installation	13
6.9	Safety precautions when operating the device	13
7	Initial operation	14
7.1	Basic equipment	14
7.1.1	TC 200 IMAGE1 S CONNECT®	14
7.1.2	TC 201 IMAGE1 S CONNECT® II	14
7.1.3	TC 300 IMAGE1 S™ H3-LINK	14
7.1.4	TC 301 IMAGE1 S™ X-LINK	14
7.1.5	TC 302 IMAGE1 S D3-LINK™	14
7.1.6	TC 304 IMAGE1 S™ 4U-LINK	14
7.2	Unpacking the equipment	14
7.3	Preparing the device for use	15
7.3.1	Installation	15
7.3.2	Connecting the potential equalization	15

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7.3.3	Connecting the IMAGE1 S™ LINK module with IMAGE1 S CONNECT®	15
7.3.4	Connecting SCB	15
7.3.5	Connecting the monitor	16
7.3.6	Connecting USB devices	16
7.3.7	Connecting the power supply	16
7.3.8	Switching on the IMAGE1 S CONNECT®	17
7.3.9	Connecting USB storage media	17
7.3.10	Connecting the camera head/video endoscope	17
7.3.11	Performing the white balance	17
7.3.12	Test for proper functioning	18
7.3.13	Installation overview (TC 200)	19
7.3.14	Installation overview (TC 201)	20
8	Operating Instructions	21
8.1	Basic information on operation	21
8.1.1	Keyboard and mouse	21
8.1.2	Camera head buttons	21
8.2	Dashboard	21
8.2.1	Description of functions of the symbols used	22
8.3	Live Menu	26
8.3.1	Description of functions of the symbols used	27
8.4	Setup Menu	32
8.4.1	General Settings (IMAGE1 S CONNECT®)	32
8.4.2	Presets	36
8.4.3	LINK-specific settings	36
8.5	Patient Data Privacy	37
9	Reprocessing	40
9.1	General safety information	40
9.2	Wipe-down disinfection of devices	41
10	Service and repair	42
10.1	Maintenance and safety check	42
10.1.1	Maintenance	42
10.1.2	Safety check	42
10.2	Servicing and repair	43
10.3	Disposal	43
10.4	Repair program	44
10.5	Limitation of liability	45
10.6	Warranty	45

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11 Technical description	46
11.1 Technical data	46
11.2 Technical documentation	49
11.3 Standard compliance	49
11.4 Directive compliance for IMAGE1 S™ (TC 200, TC 201, TC 300, TC 301, TC 302 and TC 304)	49
11.5 Software ownership and licensing	50
12 Spare parts, recommended accessories	52
12.1 Compatible camera heads for IMAGE1 S™ H3-LINK TC 300 with S technologies	52
12.2 Compatible camera heads for IMAGE1 S™ H3-LINK TC 300 without S technologies	52
12.3 Compatible camera heads for IMAGE1 S™ X-LINK TC 301 with S technologies	53
12.4 Compatible camera heads for IMAGE1 S™ X-LINK TC 301 without S technologies	53
12.5 Compatible video endoscopes for IMAGE1 S™ X-LINK TC 301 with S technologies	53
12.6 Compatible video endoscopes and camera heads for IMAGE1 S D3-LINK™ TC 302 with S technology	56
12.7 Compatible camera heads for IMAGE1 S™ 4U-LINK TC 304 with S technology	58
12.8 Recommended accessories	57
12.9 Operable SCB light sources	58
12.10 Operable SCB insufflators	58
12.11 Compatible 3D monitors	58
13 Information on electromagnetic compatibility (EMC)	59
13.1 EMC tables for IMAGE1 S™	61
14 Appendix	73
14.1 Compatible Veterinary video endoscopes for IMAGE1 S™ X-LINK TC 301 with S technologies	73
15 Subsidiaries	74

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5 General information

5.1 Description of the devices

The IMAGE1 S™ takes a unique and innovative approach to system conception. The modular design principle allows the user to create a system which is tailored to his needs. One IMAGE1 S CONNECT® or IMAGE1 S CONNECT® II module is at the heart of the IMAGE1 S™ and acts as a communication interface between the image and documentation devices.

For use, one IMAGE1 S CONNECT® or IMAGE1 S CONNECT® II plus at least one LINK module is required which receives and forwards the signal. The Light Source Control and S technologies for image optimization provide users with even more comprehensive support. An intuitive user interface assures simple and quick access as well as rapid understanding of menu navigation.

The system offers numerous state-of-the-art functions. These include:

- IMAGE1 S™ visualization modes to assist the user via contrast enhancement and/or more even brightness perception
- Light Source Control for automatic adjustment of the light source
- Integrated image and video capture for simple documentation

Functions are described in greater detail in § 8.4.

5.2 Property rights

Protected by at least one of the following US Patents: 7,821,530; 8,274,559; 7,520,853; 5,716,323; 5,788,688; 5,913,817; 6,397,286; 6,484,221; 7,212,227; US 8,089,509; US 8,199,188. For more information, visit: www.karlstorz.com.

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6 Safety instructions

**6.1 Explanation of warnings and
cautions**

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Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

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6.2 General safety information

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and all of the instruction manuals of the products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Test this device prior to each surgical procedure to ensure that it functions correctly. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

WARNING: Danger of an electric shock! This device may only be connected to a socket with protective earthing. To reduce the risk of electric shock, keep the device out of reach of patients.

WARNING: Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). The plug and cord should be routinely inspected. Do not use if damaged.

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 **WARNING:** Consult and follow the storage, transport and operating conditions prior to storing, transporting or operating, respectively. If the device is stored in a warm, humid environment, ensure that it has suitably equilibrated to the operating environment prior to use.

 **WARNING:** To reduce the risk of electrical shock, do not remove cover of unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

 **WARNING:** Perceived electrical shocks may be experienced when using the camera system in conjunction with electrosurgical units.

 **WARNING:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

 **WARNING:** Refer to the appropriate section of this manual for validated cleaning and disinfection instructions.

 **CAUTION:** Do not store liquids on or above the unit.

 **CAUTION:** In the USA, federal law restricts this device to sale or use by or on the order of a physician.

 **NOTE:** Disposal of these system products, at the end of their useful life, must be in accordance with local regulations.

 **NOTE:** Do not discard as unsorted municipal waste.

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 **NOTE:** Discard as electrical/electronic waste; recycle or reuse accordingly.

 **WARNING:** No modification of this equipment is allowed.

 **NOTE:** Consult local authorities for reuse/recycle instructions.

 **NOTE:** The presence of maximum mains voltage on SIP/SOPs is not considered a normal condition.

 **WARNING:** Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (e.g. IEC 60950 for data processing equipment and Clause 16 of 60601-1 3rd edition for medical equipment). Furthermore, all configurations shall comply with the system standard Clause 16 of 60601-1 3rd edition. Any person who connects additional equipment to the signal input part or signal output part configures a medical system, and is therefore responsible for ensuring that the system complies with the requirements of the system standard Clause 16 of 60601-1 3rd edition. If in doubt, consult the technical service department or your local representative.

 **WARNING:** When performing fluorescence imaging procedures, use suitable endoscopes and recommended light cables.

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WARNING: To ensure safe operation, do not simultaneously touch the device output connectors and the patient.

WARNING: When endoscopes are used with energized endoscopically used accessories, the patient leakage currents may be additive. This is particularly relevant when using type CF devices. In this case, a Type CF endoscopically used accessory should be used in order to minimize the total patient leakage current.

WARNING: During discharge of a defibrillator, the image displayed may take up to 4 seconds to recover and in some cases a power cycle may be required to restore the image.

WARNING: The device should be placed on a flat surface or in a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

WARNING: Failure to perform proper white balance of system may result in color distortion.

WARNING: The operator is responsible to limit physical access to the IMAGE1 S™ from unauthorized persons.

WARNING: In the case of loss of surgical image for greater than 3 seconds, a power cycle may be required to restore the image.

WARNING: When using high frequency surgical equipment, keep the working part of the active electrode in the field of view of the operator to avoid accidental high frequency burns. Avoid contact between the active electrode and metal parts of the endoscope.

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6.3 Indications for use

6.3.1 Intended use

IMAGE1 S™ is a camera control unit (CCU) for use with camera heads or video endoscopes for the visualization and documentation of endoscopic and microscopic procedures.

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6.5 Patient profile

Use of the IMAGE1 S™ is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health with regard to the endoscopic intervention and declare it adequate.

6.6 User profile (specialist physician)

- Recognized medical skills for endoscopic interventions in the relevant specialization
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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
**6.7 User profile
(assistant)**

- Completed training as medical assistant
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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6.8 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. The system is not intended for use in areas where there is a risk of explosion. This means, for example:

Wherever easily combustible and explosive inhalation anesthetics and their mixtures are used, the device may not be operated in the demarcated hazard zone. This also applies to easily combustible and potentially explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for potential equalization. It should be connected in accordance with the applicable national regulations.

6.9 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operating properly before using the device.

During treatment with the equipment, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

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7.1 Basic equipment

7.1.1 TC 200 IMAGE1 S CONNECT®

- 1 power cord, length 300 cm
- 1 DVI connecting cable, length 300 cm
- 1 SCB connecting cable, length 100 cm
- 1 3G-SDI Cable (only delivered in the USA)
- 1 KARL STORZ USB stick
- 1 Keyboard with touchpad
- 1 instruction manual

7.1.2 TC 201 IMAGE1 S CONNECT® II

- 1 power cord, length 300 cm
- 1 DVI connecting cable, length 300 cm
- 1 SCB connecting cable, length 100 cm
- 1 12G-/3G-SDI Cable
- 1 DisplayPort cable, length 300 cm
- 1 KARL STORZ USB stick
- 1 Keyboard with touchpad
- 1 instruction manual

7.1.3 TC 300 IMAGE1 S™ H3-LINK

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

7.1.4 TC 301 IMAGE1 S™ X-LINK

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

7.1.5 TC 302 IMAGE1 S D3-LINK™

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

7.1.6 TC 304 IMAGE1 S™ 4U-LINK

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

7.2 Unpacking the equipment

Carefully remove the device and accessories from the packaging. Check for missing items and evidence of damage during shipment. Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these may be useful if the unit has to be transported.

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
7.3 Preparing the device for use

7.3.1 Installation

- ① **NOTE:** Up to three IMAGE1 S™ LINK modules can be connected to the IMAGE1 S CONNECT®.
- ① **NOTE:** Do not cover the lateral air vents. Maintain a minimum distance of 6 cm to other devices.

Place the IMAGE1 S CONNECT® and the corresponding IMAGE1 S™ LINK modules on a flat surface.

7.3.2 Connecting the potential equalization

- ① **NOTE:** The unit's ground line should be installed by a qualified electrician.

The IMAGE1 S CONNECT® and the IMAGE1 S™ LINK modules are equipped with connectors for potential equalization (see fig. A and B).

1. Connect the potential equalization line to the connectors for potential equalization on all modules.
2. Connect the potential equalization line to the connector in the treatment room.

7.3.3 Connecting the IMAGE1 S™ LINK module with IMAGE1 S CONNECT®

Using the LINK cable supplied, connect the LINK module to the IMAGE1 S CONNECT® (see fig. C).

- ① **NOTE:** The IMAGE1 S CONNECT® inputs are numbered. Please connect the LINK modules in ascending order.

7.3.4 Connecting SCB

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets of the TC 200 IMAGE1 S CONNECT® or into the SCB socket of the TC 201 IMAGE1 S CONNECT® II. Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ SCB control NEO System Instruction Manual) (see fig. D).

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7.3.5 Connecting the monitor

Connect the monitor with the IMAGE1 S CONNECT® or CONNECT II (see fig. E). Depending on the types of connectors provided on the monitor, this connection can be generated using either:

- DVI-D or 3G-SDI output for IMAGE1 S CONNECT® TC 200

The primary working screen should be connected to DVI1 for TC 200.

- DisplayPort, 12G/3G-SDI or DVI-D output for IMAGE1 S CONNECT® II TC 201

The primary working screen should be connected to DP1 for TC 201.

NOTE: The connected monitor must support a resolution of 1920 x 1080 pixels.

NOTE: The connected monitor must support a resolution of 3840 x 2160 pixels to support 4K UHD.

NOTE: A compatible 3D monitor must be connected when using a D3-LINK in conjunction with 3D video endoscopes (see chapter 12.11).

7.3.6 Connecting USB devices

The USB ports of the IMAGE1 S CONNECT® are suitable for connecting up the following peripherals (see fig. F):

- Keyboard/Mouse
- Printer
- Footswitch
- USB to ACC adaptor
- USB flash drives, FAT32
- IMAGE1 PILOT TC 014

7.3.7 Connecting the power supply

Before connecting the devices to the electrical supply line, check that the values on the power sockets concur with the values on the devices' identification plates.

1. Insert the device plug of the power cord into the power supply cord receptacle of the IMAGE1 S CONNECT® (see fig. G).
2. Insert the device plug of the power cord into the power supply cord receptacle of the IMAGE1 S™ LINK modules (see fig. H).

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3. Insert the modules' power cords into the power sockets.

WARNING: Always use a hospital grade power cord with this camera.

7.3.8 Switching on the IMAGE1 S CONNECT®

Switch on IMAGE1 S CONNECT® (see fig. I).

NOTE: The connected LINK modules are turned on automatically.

7.3.9 Connecting USB storage media

Connect the USB storage media with the IMAGE1 S CONNECT® to store still images and/or videos (see fig. J).

NOTE: The USB storage medium must have a FAT 32 format.

7.3.10 Connecting the camera head/ video endoscope

Insert the camera head connector into the receptacle of the relevant LINK module (see fig. K).

NOTE: Please verify the values  for the primary and secondary sources as well as the layout beforehand (see p. 23).

NOTE: Please note the compatibility list of the respective LINK modules. (see chap. 12).

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7.3.11 Performing the white balance

NOTE: Perform a white balance before every use.

Perform the white balance by pressing the  button on the front of the IMAGE1 S CONNECT® or by pressing the  in the Live Menu.

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7.3.12 Test for proper functioning

1. Inspect **the** IMAGE1 S™ module and the camera head and/or the video endoscope for external damage.
2. Ensure that the camera head **cable** and/or the video endoscope cable is not damaged or kinked.
3. Point the camera head and/or the video endoscope towards an object and check the quality of the display on the monitor.
4. When **using** 3D video endoscopes, check to make sure the monitor's display is correct. Please ensure that **you** use the correct polarized eyewear for the monitor.

! **NOTE:** The device and accessories may only be used once the **test** for proper functioning has been successfully completed.

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


8 Operating instructions

8.1 Basic information on operation

The use of IMAGE1 S™ is divided into three menu stages:

- Dashboard (see 8.2)
- Live menu (see 8.3)
- Setup menu (see 8.4)

8.1.1 Keyboard and mouse

The keyboard and mouse can be used to navigate the menus. A keyboard is required to enter patient data.

8.1.2 Camera head buttons

All IMAGE1 FULL HD camera heads feature three buttons for retrieving programmed functions and for menu control:

- ⑮ scroll up: Scrolls up the menus, options.
- ⑰ menu/select: Displays the live menu. Once this menu is displayed, this button activates the highlighted option (identical to ENTER on a computer keyboard).
- ⑯ scroll down: Scrolls down the menus, options.

8.2 Dashboard

The Dashboard displays the devices, data storage devices and LINK modules connected to IMAGE1 S CONNECT® (see fig. L). Unconnected devices are shown in gray.

- ① **NOTE:** The Dashboard is shown when no camera head/video endoscope is connected.
- ① **NOTE:** Depending on the configuration of your system, individual icons may be deactivated or unavailable. The configuration can be modified in the Setup Menu.
- ① **NOTE:** The icons display the current settings.

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8.2.1 Description of functions of the symbols used

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


8.3 Live Menu

The live menu allows you to choose between different camera functions and change the settings of connected devices (see fig. M). Unconnected devices are shown in gray.

You can access the live menu by pressing the enter button on the keyboard or by pressing the menu button  on the camera head/video endoscope or IMAGE1 PILOT (TC014) left side button.

-  **NOTE:** The icons in the live menu are displayed as soon as a camera head is connected.
-  **NOTE:** Depending on the configuration of your system, individual icons may be deactivated or unavailable. The configuration can be modified in the Setup Menu.
-  **NOTE:** The icons display the current settings.

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8.3.1 Description of functions of the symbols used

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
NOTE: The S technologies are only available available with compatible camera heads and video endoscopes (see chapter 12).

* SPECTRA A/B is not available in the USA.

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8.4 Setup Menu

In the Setup menu, you can change the general settings of the IMAGE1 S CONNECT® as well as the link-specific settings, and optionally save them as a preset.

You can access the Setup Menu by pressing the setup button on the dashboard as well as by pressing the middle camera head button for 7 seconds in the live menu. Or press Shift-Enter on the keyboard in the setup menu.

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8.5 Patient Data Privacy

NOTE: By default the IMAGE1 S™ system does not contain any patient data. The patient data management functions (e.g. enter, access, display, store, or transmit) have to be enabled by the user.

NOTE: By enabling the patient data management functions the user accepts the responsibility for all patient related data managed with the help of IMAGE1 S™ according to applicable national data protection laws.

NOTE: KARL STORZ assumes no legal liability or responsibility for loss, correctness, damage or exposure of patient data managed with the help of IMAGE1 S™.

In the Patient Data Privacy settings, you can configure all settings related to patient information. Patient information of any kind is only available if the corresponding settings have been configured.

The configuration of the Patient Data Privacy Menu is password protected by default. The default password for the settings is: 90290245.

Pressing the Manage button opens the password management submenu (see fig. N).

Here you are offered three options (see fig. O):

- To change the current password
- To reset the password to the default password
- To reset the password protection access

Option 3 is provided to give access to the Patient Data Privacy menu settings if the current password is no longer at hand. Please note that for security reasons all patient data is removed from the device, and the device is reset to the default PDP settings.

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The Increased Patient Data Privacy item provides you with preset cyber security settings. When this setting is enabled (= active), all Patient Data Privacy sub configurations are disabled and set to the default settings. The affected sub configurations and their default values are:

▷ Display Duration

- You can determine how long the the patient data is displayed in the Patient Region of the screen.
- Default: Never

▷ Show Patient Tiles

- You can determine whether the Patient Tile will be shown on the Dashboard and Live Menu or not.
- Default: No

▷ Enable the printing of Patient Data on Prints

- You can determine whether prints will be labeled with patient information or not.
- Default: No

▷ Add Patient Names to Files Captured

- You can determine whether media files stored will be labeled with the patient information.
- Default: No

▷ Patient Data Storage

- You can specify how long data is stored on the device.
- Default: 1 day

▷ View/Save Logged Events

- You can view the Audit log by activating the Open Log button. Audit Logs can be exported via the "Save As..." function.

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The Return for Service item (see fig. P) is not affected by the Increased Patient Data Privacy item – It is always available and functions the same as the reset password protection access. **Upon confirmation, it will remove all Patient Data and reset the device settings to its default patient privacy settings.** This function is recommended for use before returning the device for service.

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9 Reprocessing

9.1 General safety information

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Do not allow water or any liquid to drip inside the unit(s).

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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9.2 Wipe-down disinfection of devices

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

① **NOTE:** KARL STORZ has qualified the resistance of the surfaces vis-à-vis wipe-down disinfection with the agent Dr. Schumacher – Cleanisept Wipes.

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10.1 Maintenance and safety check

10.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

10.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the labels relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Testing at system level in the cart includes patient leakage testing. The IEC 62353 patient leakage measurement includes testing of special (fault) conditions of mains on Applied Parts, and is conducted with one Type CF applied

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part measured in turn as per Clause 5.3.4.3.1 & Figure B2 of the IEC 62353, (limit 50 µA). When tested according to the IEC 60601-1 3.1 ed. at system level in the cart, it includes total leakage current of testing in normal conditions (limit 50 µA) as per clause # 16.6.3

Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual (see § 7.3.12).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

10.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine KARL STORZ parts only.

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Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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10.4 Repair program

For equipment, individual repair is necessary. To bridge the repair period, you usually receive a unit on loan which you then return as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ branch or authorized dealer.

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel). We reserve the right to return contaminated instruments or devices to the sender.

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- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the unit has been used in accordance with its operating instructions at all times.

10.6 Warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty specified in our Standard Terms and Conditions of Business shall additionally apply.

The device is always to be sent to your subsidiary (see chapter "Subsidiaries") even during the warranty period.

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11.2 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

11.3 Standard compliance

According to IEC 60601-1 + A1 + A2,
IEC 60601-2-18, UL 60601-1, CSA C22.2
No. 601.1-M90:

- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type CF defibrillation-proof 
- Type of protection against moisture: IPX 0

According to IEC 60601-1-2:

Please observe the information on electromagnetic compatibility in § 13 (p. 59).

11.4 Directive compliance for IMAGE1 S™ (TC 200, TC 201, TC 300, TC 301, TC 302 and TC 304)

According to Medical Device Directive (MDD):
This medical device belongs to Class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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11.5 Software ownership and licensing

With respect to KARL STORZ products containing software components, Customer has a non-exclusive, limited, non-transferable license to use the programmed logic, computer programs and/or software supplied by KARL STORZ, in connection with, and incorporated into, KARL STORZ Products System (Software) internally, only in the form in which delivered to Customer and for the sole purpose of operating in accordance with KARL STORZ's written instructions for such products sold to Customer (and for no other product). The Software, and all modifications, enhancements and upgrades thereto, will, at all times, remain the property of KSEA. Customer may not duplicate, copy, reverse-engineer, de-compile, or disassemble the Software or in any way modify the Software. Customer has no right to, and may not, create derivatives of the software, and Customer may not attempt to copy, create or re-create the source code of the Software. Any and all such modifications or enhancements to the Software by Customer, in contravention of this license, will immediately become the sole property of KSEA.

Customer hereby acknowledges and agrees (I) that the purchase of KARL STORZ products does not constitute a sale of the Software, (II) that the Software is the property of KARL STORZ, (III) that Customer neither owns nor acquires any claim or intellectual property right in or to the Software, and (IV) KARL STORZ retains all right, title, and interest in and to the Software, at all times, regardless of the form or media in or on which the original or other copies of the Software may exist. In the event of a failure of Customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the License will, without any further action by KARL STORZ or any other party, immediately terminate.

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Attribution Notice

This device contains open source software. You can view a list of the open source software included in the product, and the open source software license terms, on the IMAGE1 S CONNECT® device itself. From the Dashboard, select "Setup/General", then "Software Licenses".

To obtain a copy of the source code being made publicly available by KARL STORZ Imaging, Inc. ("KARL STORZ") related to open source software used in this KARL STORZ product, you should send your request, along with product part number and serial number, in writing by email to opensource@karlstorz.com or by postal mail to:

KARL STORZ Imaging, Inc.
Attention: Open Source Software Inquiries
1 S. Los Cameros
Goleta, CA 93117
USA

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**Spare parts,
recommended accessories**

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**12 Spare parts,
recommended accessories**

**12.1 Compatible camera heads
for IMAGE1 S™ H3-LINK
TC 300 with S technologies**

Name	Article number
IMAGE1 S™ H3-Z Three- Chip FULL HD Camera Head	TH 100
IMAGE1 S™ H3-ZI Three- Chip FULL HD Camera Head	TH 101
IMAGE1 S™ H3-Z FI Three- Chip FULL HD Camera Head	TH 102
IMAGE1 S™ H3-P Three- Chip FULL HD Pendulum Camera Head	TH 103
IMAGE1 S™ H3-ZA Three- Chip FULL HD Camera Head	TH 104
IMAGE1 S™ H3-M COMIEW® Three-Chip FULL HD Camera Head	TH 106

**12.2 Compatible camera heads
for IMAGE1 S™ H3-LINK
TC 300 without S technologies**

Name	Article number
IMAGE1 H3-Z Three-Chip FULL HD Camera Head	22 2200 55-3
IMAGE1 H3-ZI Three-Chip FULL HD Inline Camera Head	22 2200 53-3
IMAGE1 H3-Z FI Three-Chip ICG FULL HD Camera Head	22 2200 85-3*
IMAGE1 H3-P Three-Chip FULL HD Pendulum Camera Head	22 2200 56-3
IMAGE1 H3-ZA Three-Chip FULL HD Camera Head	22 2200 61-3
IMAGE1 H3-FA Three-Chip FULL HD Camera Head	22 2200 60-3

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WARNING: IMAGE1 S™ should not be placed in the immediate vicinity of other electrical equipment. If this cannot be avoided, please ensure that the performance of either product is not affected by unintended electromagnetic coupling.

WARNING: IMAGE1 S™ is only intended for use by physicians and medical assistants. This is a CISPR Class A medical equipment/system. In a domestic environment, IMAGE1 S™ may cause radio interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the unit or filtering the connection to the public mains network.

NOTE: The EMC tables and other guidelines that are included in this Instruction Manual provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

NOTE: This device has been tested and complies with the EMC limits values for the Medical Device Directive 93/42/EEC [EN 55011 Class A and IEC 60601-1-2].

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These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The equipment generates uses and can radiate radio frequency energy. If not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference with other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the individual devices
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected
- Consult the manufacturer or field service technician for help.

WARNING: The use of portable and mobile RF equipment may have impact on this and other pieces of medical equipment.

WARNING: The use of an accessory or cable with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the IMAGE1 S™. When using an accessory or cable with IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

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System Cables and Maximum Lengths used for EMC Compliance				
Type	Shielded	Max Length [m]	Ferrite	Used for
LINK cable	Yes	2	No	Interconnection of inter-module LINK signals
BNC to BNC	Yes	10	No	Connection of 3G/12G-SDI video signals
DVI	Yes	3	No	Connection of DVI signals
DisplayPort (DP)	Yes	3	No	Connection of DisplayPort signals
4-pin mini-DIN	Yes	1	No	Connection of SCB signals
USB	Yes	2	No	Connection to printer, storage device, keyboard
Power cord	No	3	No	Connection of Device to AC Mains
Ethernet	No	2	No	Connection to LAN (RJ-45) port
IMAGE1 PILOT TC 014	Yes	5	No	Connection to VITOM® 3D
External synchronization cable	No	2.5	No	Connection to compatible light source

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Table 1		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic enforcement – guidance
RF emissions acc. to CISPR 11	Group 1	The IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The user of IMAGE1 S™ should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Air discharge	± 6 kV Contact discharge ± 8 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Rapid transient electrical interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	± 2 kV/1 kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
* Note: U_T is the mains alternating voltage before applying the test levels.			

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Table 4 Guidance and manufacturer's declaration – electromagnetic immunity			
The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ should ensure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons. ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IMAGE1 S™ is used exceeds the applicable RF compliance level above, the IMAGE1 S™ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the IMAGE1 S™. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 6

Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™

The IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

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14 Appendix

14.1 Compatible Veterinary video
endoscopes for IMAGE1 S™
X-LINK TC 301
with S technologies

Name	Article number (PAL/NTSC)
Specialty VET	
Veterinary Video Endoscope	60914 NKS/PKS
Veterinary Video Endoscope	60714 NKS/PKS
Veterinary Video Endoscope	60511 NKS/PKS
Veterinary Video Endoscope	60332 NKS/PKS
Veterinary Video Endoscope	60130 NKS/PKS
Veterinary Video Endoscope	60125 NKS/PKS
Veterinary Video Endoscope	60118 NKS/PKS

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@ksesa.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N.W. 57th Avenue, Suite 490
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No. 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 Mexico D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing America Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar — Raim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufrautegui 627 6° Piso, 81639 CAA — Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østlensveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kunoors Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Tavaihe 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68587000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich — Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-it@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Caprašva 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gonice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 8113-1188, Fax: +86 21 8113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610044, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karsltoz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karsltozvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S A de C.V.
Lago Constanza No 326, Col. Chapulterec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America do Sul Ltda.
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, Sao Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karsltoz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufnategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karsltoz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvi 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsatragrand 14
12724 Skarholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karsltoz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4. Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 96824775
E-Mail: asiakaspalvelu@karsltoz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karsltoz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karsltoz-uk.com

KARL STORZ Endoscopie Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopleer.be

KARL STORZ Endoscopie France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karsltoz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: stolz-austria@karsltoz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karsltoz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell' Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karsltoz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karsltoz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) – Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai – United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karsltoz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 – Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karsltoz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karsltoz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karsltoz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karsltoz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karsltoz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karsltoz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karsltoz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karsltoz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karsltoz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karsltoz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karsltoz.co.jp

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Important note for users of KARL STORZ devices and instruments

1 Important note for users of KARL STORZ devices and instruments

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Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual, in combination with the IMAGE1 S™ instruction manual (96206286DER), is intended to aid you in the proper installation, connection, and operation of this device. All required details and all actions on your part are clearly presented and explained. Please read this manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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1 Important note for users of KARL STORZ devices and instruments	IV
2 Safety instructions	3
2.1 Explanation of warnings and cautions	3
2.2 General safety information	4
3 Intended Use Camera Heads	6
3.1 Contraindications	6
3.2 User qualification	6
4 Images of the equipment	7
5 Symbols employed	8
5.1 Symbols for operation	8
5.2 Label and packaging symbols	8
6 Operating instructions for the IMAGE1 S™ HX FI camera heads	9
6.1 Frequency of the IMAGE1 S™ HX FI camera heads (50/60 Hz)	9
6.2 Operating the IMAGE1 S™ HX FI camera heads	9
6.2.1 Connecting an endoscope and light cable	9
6.2.2 Adjusting the focus	10
6.3 Components of the KARL STORZ FI System	10
7 Technical description	11
7.1 Specifications	11
7.2 Storage and Operating conditions	11
7.3 Standard compliance	12
7.4 Directive compliance	12
8 Reprocessing (not validated for US customers)	13
8.1 Care and handling	13
8.2 Reprocessing	14
8.2.1 Reference chart	14
8.2.2 General warnings	15
8.3 Accessories	17
8.4 Preparation for cleaning and disinfection	17
8.4.1 Manual wipe-down disinfection	17
8.4.2 Manual precleaning	18
8.4.2.1 Brushing the surfaces	18
8.4.2.2 Auxiliary pretreatment with ultrasound	18

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8.4.3	Manual cleaning.....	18
8.4.4	Manual disinfection.....	18
8.5	Assembly, inspection and care.....	19
8.6	Packaging systems.....	20
8.7	Sterilization.....	21
8.7.1	Hydrogen peroxide (H ₂ O ₂)- sterilization – ASP STERRAD®.....	21
8.7.2	Hydrogen peroxide (H ₂ O ₂)- sterilization – STERIS® AMSCO® V-PRO™1.....	22
8.7.3	Ethylene oxide sterilization (EO).....	22
8.7.4	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®.....	23
8.7.5	High-Level Disinfection.....	24
8.7.6	CIDEX (2.4 % glutaraldehyde solution).....	24
8.7.7	Resert XL HLD (2.0 % hydrogen peroxide solution).....	25
8.8	Reprocessing limits.....	26
9	Maintenance.....	27
9.1	Maintenance and safety check.....	27
9.1.1	Maintenance.....	27
9.1.2	Safety check.....	27
9.2	Servicing and repair.....	28
9.3	Disposal.....	28
9.4	Repair program.....	29
9.5	Limitation of liability.....	30
9.6	Warranty.....	30
10	Compatible camera control unit (CCU).....	31
11	Spare parts, recommended accessories.....	32
12	Information on electromagnetic compatibility (EMC).....	33
13	Subsidiaries.....	45

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2.2 General safety information

 **WARNING:** Read this instruction manual thoroughly before using the equipment.

 **WARNING:** Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.

 **WARNING:** Perceived electric shocks may be experienced when using the camera system in conjunction with electrosurgical units.

 **WARNING:** When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.

 **WARNING:** When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. It should be noted that this is particularly important if a Type CF applied part endoscope is used, in which case a Type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

 **WARNING:** Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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WARNING: All endoscopic devices (e.g., HF surgical probes) used in conjunction with KARL STORZ endoscopic cameras must comply with IEC 60601-2-18, the standard for medical electrical equipment – particular requirements for the safety of endoscopic equipment. The insulation of the endoscopic devices must be designed to withstand the highest HF voltage encountered during surgery.

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WARNING: Test all equipment prior to each surgical use. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

WARNING: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

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3 Intended Use Camera Heads

The camera head is used together with the camera control unit (CCU) and is attached to a telescope during endoscopic or microscopic procedures. In addition, the camera head is used for Fluorescence Imaging applications.

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3.2 User qualification

The IMAGE1 S™ HX FI camera heads may only be used by physicians and medical assistants who have a specialized qualification and who have been instructed in use of the IMAGE1 S™ HX FI camera heads.

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

- ① Instrument coupler
- ② Focusing ring/knob
- ③ User control button: for scrolling (downwards) through the menu or activating a camera function
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑤ User control button: for accessing the menu/ making a selection

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6 Operating instructions for the IMAGE1 S™ HX FI camera heads

6.1 Frequency of the IMAGE1 S™ HX FI camera heads (50/60 Hz)

The IMAGE1 S™ HX FI camera heads are suitable for operation at either 50 Hz or 60 Hz. The camera head must be selected to suit the local conditions.

! **NOTE:** The method for switching frequency is described in the instruction manual IMAGE1 S™ – TC 200, TC 300, TC 301 (96206286DER).

6.2 Operating the IMAGE1 S™ HX FI camera heads

6.2.1 Connecting an endoscope and light cable

The IMAGE1 S™ HX FI camera heads have an integrated instrument coupler and integrated telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the coupler clockwise and insert the endoscope eyepiece. When connecting the endoscope to a Pendulum camera head, the Pendulum lock must click into place before the outer ring of the grasping mechanism is turned. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the endoscope eyepiece by tightening the knurled screw on the light cable through a quarter turn.

! **CAUTION:** Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

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- ① **NOTE:** For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (96206342ZS), which can be downloaded or requested at www.karlstorz.com.
- ② **NOTE:** For U.S. Customers, please refer to PI-0000011 'Fiber Optic Light Cables Care and Maintenance'.

6.2.2 Adjusting the focus

The focus on HX FI camera heads can be adjusted by turning the focusing ring/knob ② on the camera lens (see Fig. 1 and 2).

6.3 Components of the KARL STORZ FI System

The KARL STORZ FI System comprises the KARL STORZ D-LIGHT C/AF/P system including suitable KARL STORZ fluid light cables, KARL STORZ FI endoscope or FI videoscope and KARL STORZ FI camera with coordinated KARL STORZ camera control unit.

Refer to the system description for KARL STORZ Fluorescence Imaging (FI) Systems (see instruction manual 96206144DER) for the recommended configuration.

- ① **NOTE:** Please observe the instruction manuals for the KARL STORZ camera control units IMAGE1 S™ (cat. no. 96206286DER).

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**8 Reprocessing (not validated
for US customers)**

In the following chapter, you will find information on care and cleaning of the IMAGE1 S™ HX FI camera heads.

8.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable "kinks" may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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8.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile upon delivery. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, national regulations regarding personal safety must be observed.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ.

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NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

CAUTION: Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be reprocessed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

CAUTION: Chemical residue on the camera head plug can cause failure of the camera electronics.

CAUTION: The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.

CAUTION: The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

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8.4.2 Manual precleaning

8.4.2.1 Brushing the surfaces

Depending on the degree of soiling and the degree to which residues have dried, heavy soiling must be cleaned from the surfaces of instruments with the aid of a brush/sponge under cold running water until any visible contamination has **been** completely removed.

8.4.2.2 Auxiliary pretreatment with ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

8.4.3 Manual cleaning

The medical device must be completely immersed in a cleaning solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a brush or sponge. The device must then be rinsed with cold water to ensure neutralization.

8.4.4 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces and joints are dried completely with (preferably sterile) compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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The following sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ HX FI camera heads:

8.7.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of dimensions and material.

CAUTION: For more detailed information, please consult the user handbook of the respective device.

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NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S
- STERRAD® NX® 'standard' cycle
- STERRAD® 100NX® 'standard' cycle
- 'DUO' cycle

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**8.7.2 Hydrogen peroxide (H₂O₂)
sterilization – STERIS® AMSCO®
V-PRO™1**

Detailed information on the selection of the appropriate cycle for the various device generations is available from the STERIS® manufacturer.

! **CAUTION:** V-PRO™ sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

8.7.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100 % ethylene oxide at 55 °C and a hold time of 30 – 45 minutes.

! **WARNING:** Medical devices must be dismantled for sterilization.

i **NOTE:** Place the camera head in a sterilization tray. Do not place any instruments on top of the camera.

! **CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing technique of the EO sterilizer (concentration, process control).

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

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**Reprocessing
(not validated for US
customers)**

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To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head with sterile water is **very** important for removing any residual disinfectant.

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9.1 Maintenance and safety check

9.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

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Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.

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Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual (see Chapter 7.3.12 of the IMAGE1 S™ CCU instruction manual no. 96206286DER).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

9.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

9.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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9.5 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

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Item	Cat. No.
IMAGE1 S CONNECT™ in combination with IMAGE1 S™ X-LINK	TC 200 TC 301

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Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Plastic container for sterilization and storage of IMAGE1 S™ HX FI (TH 112) and IMAGE1 S™ HX-P FI (TH 113) camera heads	39301 HXTS
Instruction Manual	96286002DER
Instruction manual IMAGE1 S™	96206286DER

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12 Information on electromagnetic compatibility (EMC)

The CGU and camera heads were tested as a system.

WARNING: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Please observe the EMC instructions in this instruction manual during installation and operation.

WARNING: HX camera heads are only intended for use by physicians and medical assistants. This is a CISPR Class A medical equipment/system. In a domestic environment, HX camera heads may cause radio interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the unit or filtering the connection to the public mains network.

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WARNING: The use of an accessory or cable with the HX camera heads other than those specified in this manual may result in increased emissions or decreased immunity of the HX camera heads.

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Table 1		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The HX camera heads are intended for use in the electromagnetic environment specified below. The customer or user of the HX camera heads should ensure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic enforcement – guidance
RF emissions acc. to CISPR 11	Group 1	The HX camera heads use RF energy only for their internal function. The customer or user of the HX camera heads should make sure that they are used in such an environment. Therefore, their RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The HX camera heads are suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The HX camera heads are intended for use in the electromagnetic environment specified below. The user of HX camera heads should make sure that they are used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Air discharge	± 6 kV Contact discharge ± 8 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Rapid transient electrical interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	± 2 kV/1 kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV Differential mode ± 2 kV Common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HX camera heads requires continued operation during power mains interruptions, it is recommended that the HX camera heads be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the HX camera heads further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
 *NOTE: U_T is the mains alternating voltage before applying the test levels.			

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Table 4			
Guidance and manufacturer's declaration – electromagnetic immunity			
The HX camera heads are intended for use in the electromagnetic environment specified below. The customer or user of the HX camera heads should ensure that they are used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the HX camera heads, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range [°]. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.			
^a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. [°] The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges. [°] Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HX camera heads are used exceeds the applicable RF compliance level above, the HX camera heads should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the HX camera heads. ^a Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 6

**Recommended separation distances between portable and mobile RF communications equipment and
the HX camera heads**

The HX camera heads are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the HX camera heads can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the HX camera heads as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.6 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4190 (Canada only)
E-Mail: info.canada@karl-storz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseaa.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cerritos Road
Galeita, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karl-storzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia Mexico S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No. 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karl-storz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karl-storz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karl-storz.co.za

TOO KARL STORZ Endoscopy Kazakhstan
6, Saryarka str. 8C, Arman, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karl-storz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A - Unit 7WA - 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karl-storz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karl-storz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karl-storz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karl-storz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karl-storz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jiangguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 56381888, Fax: +86 10 56381999
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linxiang Road, Wuhou District,
Chengdu, 610044, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyang District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1261, Fax: +86 20 8732-1286
E-Mail: info@karl-storz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karl-storz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karl-storz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
238 Jeonggil-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karl-storz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 9672 6399
E-Mail: mingwang@karl-storz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 31745-00
E-Mail: philippines@karl-storz.com

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karl-storz.co.jp

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseacom

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 53th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 7204 1097, Fax: +53 7204 1098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapultec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America do Sul Ltda.
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, São Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsensvei 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsåstragård 14
12724 Skårholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4. Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stonler.be

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) - Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-I
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 - Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jelion Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: ningwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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**Important note for users of
KARL STORZ devices and
instruments**

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**1 Important note for users of
KARL STORZ devices and
instruments**

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

1.1 Intended use

The camera head is used together with the camera control unit and is connected to a telescope during endoscopic or microscopic interventions.

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1 Important note for users of KARL STORZ devices and instruments	III
1.1 Intended use	III
2 Safety instructions	3
2.1 Explanation of warnings and cautions	3
2.2 General safety information	4
3 Images of the equipment	6
4 Symbols employed	7
4.1 Symbols for operation	7
4.2 Label and packaging symbols	7
5 Operating instructions	8
5.1 Frequency of the High Definition camera head (50/60 Hz)	8
5.2 Using endoscopy camera heads (TH100, TH101, TH103, TH104)	8
5.2.1 Connecting an endoscope and light cable	8
5.2.2 Adjusting the focus	9
5.2.3 Adjusting the Image size with zoom camera heads	9
5.2.4 Operating the user control buttons	10
5.3 Using the microscope camera head TH106 IMAGE1 S™ H3-M COVIEW®	11
5.3.1 Detaching the camera head cable	11
5.3.2 Connecting the camera head to the TV adaptor and surgical microscope	11
5.3.3 Connecting the camera cable	12
5.3.4 Operating the user control buttons	12
6 Technical description	13
6.1 Specifications	13
6.2 Storage and operating conditions	14
6.3 Standard compliance	14
6.4 Directive compliance	14
6.5 Information on electromagnetic compatibility (EMC)	14
7 Reprocessing	15
7.1 Care and handling	15
7.2 Reprocessing	16
7.2.1 Reference chart	16
7.2.2 General warnings	17
7.3 Preparation for cleaning and disinfection	19
7.3.1 Manual wipe-down disinfection	19
7.3.2 Manual cleaning	20
7.3.3 Manual disinfection	20

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7.4	Inspection and care.....	21
7.5	Packaging systems.....	21
7.6	Sterilization.....	21
7.6.1	Steam sterilization in the fractionated prevacuum procedure only suitable for autoclavable camera heads.....	22
7.6.2	Steam sterilization using the gravitation method only suitable for autoclavable camera heads.....	22
7.6.3	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®.....	23
7.6.4	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO®.....	24
7.6.5	Ethylene oxide sterilization (EO).....	25
7.6.6	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1 (except U. S.).....	25
7.6.7	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E (only valid in U. S.).....	25
7.6.8	High-Level Disinfection.....	25
7.6.9	CIDEX.....	26
7.6.10	Resert XL HLD.....	26
7.7	Reprocessing limit.....	27
8	Maintenance.....	28
8.1	Maintenance and safety check.....	28
8.1.1	Maintenance.....	28
8.1.2	Safety check.....	28
8.2	Servicing and repair.....	29
8.3	Disposal.....	29
9	Repairs and warranty.....	30
9.1	Repair program.....	30
9.2	Limitation of liability.....	31
9.3	Warranty.....	31
10	Compatible camera control units.....	32
11	Spare parts, recommended accessories.....	33
12	Subsidiaries.....	34

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2 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

2.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

This instruction manual is intended to serve only as a supplement to the KARL STORZ video systems instruction manual, which was enclosed with your camera control unit. To ensure the safe and efficient use of this KARL STORZ product, it is essential that you read the entire video system instruction manual before connecting the camera system.

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WARNING: All endoscopic devices (e.g., HF surgical probes) used in conjunction with KARL STORZ endoscopic cameras must comply with IEC/EN 60601-2-18, the standard for medical electrical equipment – particular requirements for the safety of endoscopic equipment. The insulation of the endoscopic devices must be designed to withstand the highest HF voltage encountered during surgery.

HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.

WARNING: Risk of blindness. Never look into the open end of a light cable.

WARNING: Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.

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



- ① Instrument coupling,
C-MOUNT connection with H3-M
- ② Focusing ring/knob
- ③ Zoom ring
- ④ User control button: for scrolling (downwards)
through the menu or activating a camera function
- ⑤ User control button: for scrolling (upwards)
through the menu or activating a camera function
- ⑥ User control button: for accessing the menu/
making a selection

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4 Symbols employed

4.1 Symbols for operation

	Type CF device/Defib. resistant
	CE marking
	Zoom in
	Zoom out

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4.2 Label and packaging symbols

	Fragile, handle with care
	Serial number
	Number of products in the product packaging
	CE marking
	This side up
	Consult instructions for use
	Manufacturer

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5 Operating instructions

5.1 Frequency of the High Definition camera head (50/60 Hz)

The IMAGE1 S™ H3 camera heads are suitable for use with frequencies of both 50 Hz and 60 Hz. The frequency can be adapted to suit local conditions.

NOTE: The method for switching frequency is described in the instruction manual IMAGE1 S™ – TC 200, TC 300, TC 301 (96206286DER).

5.2 Using endoscopy camera heads (TH 100, TH 101, TH 103, TH 104)

5.2.1 Connecting an endoscope and light cable

IMAGE1 S™ H3 camera heads have an integrated instrument coupler and telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the instrument coupler clockwise and insert the eyepiece of the endoscope into the coupler. When connecting the endoscope to a Pendulum camera head, the Pendulum lock must click into place before the outer ring of the grasping mechanism is turned. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the eyepiece of the endoscope by tightening the knurled screw on the light cable through a quarter turn.

CAUTION: Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

NOTE: For further information regarding light cables, please refer to the manual "Fiberoptic light cables" (97000213), which can be downloaded or requested at www.karlstorz.com.

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5.2.2 Adjusting the focus

The focus can be adjusted by turning the focusing ring/knob on the camera lens.

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
5.2.3 Adjusting the image size with zoom camera heads

The image size can be adjusted by turning the zoom ring.

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5.2.4 Operating the user control buttons

- ④ Scrolling downwards: Scroll downwards through menus, options, etc. Can also be assigned a programmable function/feature.
- ⑤ Scrolling upwards: Scroll upwards through menus, options, etc. Can also be assigned a programmable function/feature.
- ⑥ Menu/select: Access the menu. Once the menu is displayed, this button can be used to select the highlighted option (identical to Enter on a computer keyboard). Hold this button down for seven seconds to call up the "Options" menu.

For detailed instructions on the programming and use of the camera head buttons, please refer to the KARL STORZ IMAGE1 S™ (Art. No.: 96206286DER) and KARL STORZ IMAGE1 (Art. No.: 96206114DER) instruction manuals.

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5.3 Using the microscope camera head TH 106 IMAGE1 S™ H3-M COVIEW®

The IMAGE1 S™ H3-M COVIEW® microscope camera head is connected to a suitable interface on the surgical microscope using a TV adaptor. Suitable interfaces are optical beam splitters and co-observer connectors. These options are available for all surgical microscopes.

Some surgical microscopes are equipped with a C-MOUNT connection that is directly coupled with the IMAGE1 S™ H3-M COVIEW® camera head.

5.3.1 Detaching the camera head cable

The IMAGE1 S™ H3-M COVIEW® camera head is delivered with a 9 m cable as standard. Detaching the camera cable from the camera head is recommended in order to make installation easier. To detach the camera cable, unscrew the cable connector from the camera head by turning in the counterclockwise direction and carefully pull out the cable.

5.3.2 Connecting the camera head to the TV adaptor and surgical microscope

To connect the the IMAGE1 S™ H3-M COVIEW® camera head to the surgical microscope, a TV adaptor is required (e.g. KARL STORZ QUINTUS® TV adaptor). Screw the camera head carefully onto the C-MOUNT connection of the TV adaptor and then secure the camera head and TV adaptor to the co-observer interface of the surgical microscope. For more detailed installation instructions, please refer to the instruction leaflet provided with the QUINTUS® TV adaptor.

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

5.3.3 Connecting the camera cable

To connect the camera head cable to the camera head, align the tip of the cable plug carefully with the recess provided and press the cable carefully onto the camera head. Screw the cable connector onto the camera head by turning in the clockwise direction.

CAUTION: To avoid damage, the camera cable should not be connected to the Camera Control Unit (CCU) during installation!

CAUTION: The camera cable should NOT be connected to/disconnected from the camera head when the CAMERA CONTROL UNIT (CCU) is switched on.

5.3.4 Operating the user control buttons

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For detailed instructions on the programming and operation of the camera head buttons, please refer to the KARL STORZ IMAGE1 instruction manual.

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6.2 Storage and operating conditions

Refer to the compatible camera control unit (CCU) manual for relevant storage and operating conditions. The CCU and camera head were evaluated as a system for IEC 60601-1.

6.3 Standard compliance

According to IEC 60601-1; IEC 60601-2-18; UL 60601-1

- Protection of applied part against electric shock: Type CF 

6.4 Directive compliance

According to the Medical Device Directive (MDD):

This medical device belongs to Class I

This medical device bears a CE mark in accordance with MDD 93/42/EEC.

6.5 Information on electromagnetic compatibility (EMC)

Refer to the Electromagnetic Compatibility (EMC) information Appendix of the compatible camera control unit (CCU) manual for relevant electromagnetic compatibility (EMC) information.

The EMC warning statements, caution, notes and Emission/Immunity limits specified in the EMC Tables of the compatible CCU apply to the camera head. The CCU and camera head were tested as a system.

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7.2.2 General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

 **WARNING:** Risk of infection: These medical devices are not sterile upon delivery. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

 **WARNING:** When carrying out any work on contaminated medical devices, national regulations regarding personal safety must be observed.

 **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

 **CAUTION:** Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ.

 **CAUTION:** National laws and regulations must be observed.

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i NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

! CAUTION: Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be processed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

! CAUTION: Chemical residue on the camera head plug can cause failure of the camera electronics.

! CAUTION: The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.

! CAUTION: The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

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7.3 Preparation for cleaning and disinfection

Immerse the camera head in a suitable container containing a pH-neutral (pH value of 6.0 – 8.0) enzymatic cleaning solution as per the manufacturer's instructions immediately after use to prevent blood, protein and other contaminants from drying onto the equipment.

CAUTION: Inspect the camera head cable for breaks and cuts. Camera heads with damaged cables must not be immersed. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

Heavy soiling must be removed from the camera head immediately after use. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

7.3.1 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material incompatibility must be observed.

CAUTION: The TH 106 IMAGE1 S™ H3-M COVIEW® microscope camera head may only be cleaned by wipe-down disinfection. Other reprocessing or sterilization methods are not permitted.

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7.3.2 Manual cleaning

The camera head and camera head cable must be completely immersed, preferably in a pH-neutral enzymatic cleaning solution as per the manufacturer's instructions and microbiologically pure/completely demineralized water. It must be ensured that the surfaces are covered and that no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with a soft brush (Art. No. 27652) or sponge. Finally, it must be rinsed with microbiologically pure, sterile or completely demineralized water to ensure neutralization. Thorough rinsing of the camera head is necessary for removing any soiling and cleaning agent residues which could interfere with disinfection.

7.3.3 Manual disinfection

The camera head and the camera head cable must be completely immersed in a disinfectant solution. It must be ensured that the surfaces are covered and that no air bubbles are present. Following exposure, rinse the camera head three times with microbiologically pure water, sterile water, or completely demineralized water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Thorough rinsing of the camera head is necessary for removing any residual disinfectant that could interfere with sterilization.

All the surfaces are then dried completely with sterile compressed air. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose. Wipe the camera head connector with 70% isopropyl alcohol to remove any residual disinfectant.

Do not allow exposed glass windows to air dry. 70% isopropyl alcohol may be applied to glass surfaces with a cotton tip applicator to prevent streaks and spots. Dry the surfaces thoroughly with a cotton tip applicator after applying the alcohol.

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7.4 Inspection and care

The cleaned camera head and cleaned camera head cable must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the camera head and the camera head cable must be manually cleaned and subjected to a full cleaning procedure once more.
- Camera heads with damaged cables must be withdrawn from use.

CAUTION: Inspect the camera head cable for breaks and cuts. Camera heads with damaged cables must not be sterilized or disinfected. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

- Afterwards, a functional check must be carried out.

Before sterilization and/or disinfection, coil the camera head cable into loops with a diameter of at least 15 cm (6"). In doing so, do not kink or twist the cable.

7.5 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 - 10, EN ISO 11607 Parts 1 + 2, DIN 58953). National regulations must be observed.

7.6 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments". The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

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The following sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ H3 camera heads:

7.6.1 Steam sterilization in the fractionated prevacuum procedure only suitable for autoclavable camera heads

! **CAUTION:** Of the camera heads compatible with IMAGE 1 HUB™ HD and IMAGE1 S™ only, IMAGE1 S™ H3-ZA (TH 104) is autoclavable.

i **NOTE:** Place the camera head in a suitable sterilization container (e.g. 39301 Z3TS).

The autoclavable camera head and the camera head cable must be sterilized using a fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134 °C – 137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

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7.6.2 Steam sterilization using the gravitation method only suitable for autoclavable camera heads

! **CAUTION:** Of the camera heads compatible with IMAGE 1 HUB™ HD and IMAGE1 S™ only, IMAGE1 S™ H3-ZA (TH 104) is autoclavable.

i **NOTE:** Place the camera head in a suitable sterilization container (e.g. 39301 Z3TS).

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The autoclavable camera head and the camera head cable must be sterilized using a gravitation procedure at 121 °C – 122 °C with an exposure time of 30 minutes.

It is recommended to leave the camera heads to dry for at least 4 minutes, so that the card-edge connector is completely dry before its insertion.* After steam sterilization, the camera heads should not be touched until adequately cooled.

! **CAUTION:** Do not immerse the camera head in cold liquid after removing it from the steam sterilizer. Due to the "temperature shock", this may result in damage to the camera lens and invalidation of the warranty.

* When using steam sterilizers without an integrated drying cycle, the card-edge connector can be dried off with a lint-free, sterile cloth once it has cooled down. Before insertion, all excess moisture must be removed from the card-edge connector.

7.6.3 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

! **WARNING:** Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of dimensions and material.

! **CAUTION:** For more detailed information, please consult the user handbook of the respective device.

i **NOTE:** To ensure that the medical device concerned can be sterilized using the different STERRAD® units, refer to the "STERRAD® Sterility Guide" (www.sterradsterilityguide.com).

! **CAUTION:** STERRAD® sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

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The following STERRAD® sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ H3 camera heads:

- STERRAD® 100S
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
„short“ cycle
- TH 103 H3-P
„long“ cycle
- STERRAD® NX
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
„standard“ cycle
- STERRAD® 100NX
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 S H3-ZA
„standard“ cycle
„DUO“ cycle

7.6.4 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO®

More detailed information regarding the cycle used for various device generations can be obtained from the manufacturer STERIS®.

CAUTION: V-PRO™ sterilization may cause cosmetic changes to the camera head that do not usually impact the camera's functionality.

The following STERIS® AMSCO® sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ H3 camera heads:

- V-PRO™ 1
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
Equivalent to „lumen“ Cycle
- V-PRO™ 1 Plus
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
„lumen“ cycle
- V-PRO™ maX
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
„lumen“ cycle

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7.6.5 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100 % ethylene oxide at 55 °C and a hold time of 30-45 minutes.

 **NOTE:** Place the camera head in a sterilization tray. Do not place any instruments on top of the camera.

 **CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated for the materials by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing procedure of the EO sterilizer (concentration, process control).

7.6.6 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1 (except U. S.)

More detailed information regarding sterilization parameters can be obtained from the manufacturer STERIS®.

7.6.7 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E (only valid in U. S.)

More detailed information regarding sterilization parameters can be obtained from the manufacturer STERIS®.

7.6.8 High-Level Disinfection

 **WARNING:** High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

 **CAUTION:** Any deviations from the recommended parameters for disinfection must be validated by the user.

 **CAUTION:** Before disinfection, the instruments must be thoroughly cleaned, rinsed, and dried.

 **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life.

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The following high-level disinfection methods have been validated and approved by KARL STORZ for the IMAGE1 S™ H3 camera heads:

7.6.9 CIDEX (2.4% glutaraldehyde solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.4% glutaraldehyde solution for a minimum of 45 minutes at 25 °C (77°F). It must be ensured that the surfaces are covered and that no air bubbles are present.

! CAUTION: Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head and the card-edge connector with sterile water is very important for removing any residual disinfectant.

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7.6.10 Resert XL HLD (2.0% hydrogen peroxide solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20 °C (68 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

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9 Repairs and warranty

9.1 Repair program

CAUTION: Opening the camera head will result in damage and void the warranty.

To bridge the repair period, you will, as a rule, receive a unit on loan which must then be returned to KARL STORZ as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel).

We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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9.2 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

9.3 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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NOTE: The S technologies are available with IMAGE1 S™.

* IMAGE1 camera control units are available with various configurations and image outputs. These camera control units are compatible with IMAGE1 S™ camera heads.

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**11 Spare parts,
recommended accessories**

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Extension cable 5.2 m	22 2200 72
Plastic Container for Sterilization and Storage of camera heads IMAGE1 S™ H3-Z (TH 100) and IMAGE1 S™ H3-ZA (TH 104), autoclavable	39301 Z3TS
H3-M COVIEW® Camera Head Cable 9 m	22 2200 71
Instruction Manual	96206511DER
Plastic Container for Sterilization and Storage of camera heads IMAGE1 S™ H3-ZI (TH 101) and IMAGE1 S™ H3-P (TH 103)	39301 PHTS

* For safe and comfortable working with KARL STORZ cameras we recommend the use of sterile drapes.

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kazakhstan
6, Saryarka str. BC "Arman", office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A - Unit 7WA - 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
11th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25/F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18/F Beijing IFC
No. 8, Jiangguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguo Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Lingang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
J-Park Bundang1, 102dong, 8F
239 Jeongail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mmgwan@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Ujima
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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1 Important note for users of KARL STORZ devices and instruments.....	II
2 Safety instructions.....	4
2.1 Explanation of warnings and cautions.....	4
2.2 General safety information.....	5
3 Intended Use Camera Heads.....	7
3.1 Contraindications.....	7
3.2 User qualification.....	7
4 Images of the equipment.....	8
5 Symbols employed.....	9
5.1 Symbols for operation.....	9
5.2 Label and packaging symbols.....	9
6 Operating instructions for the IMAGE1 S™ 4U camera head.....	10
6.1 Frequency of the IMAGE1 S™ 4U camera head (50/60 Hz).....	10
6.2 Operating the IMAGE1 S™ 4U camera head.....	10
6.2.1 Connecting an endoscope and light cable.....	10
6.2.2 Adjusting the focus.....	11
7 Technical description.....	11
7.1 Specifications.....	11
7.2 Storage and Operating conditions.....	11
7.3 Standard compliance.....	12
7.4 Directive compliance.....	12
8 Reprocessing (not validated for US customers).....	13
8.1 Care and handling.....	13
8.2 Reprocessing.....	14
8.2.1 Reference chart.....	14
8.2.2 General warnings.....	15
8.3 Accessories.....	17
8.4 Preparation for cleaning and disinfection.....	17
8.4.1 Manual wipe-down disinfection.....	17
8.4.2 Manual precleaning.....	18
8.4.2.1 Brushing the surfaces.....	18
8.4.2.2 Auxiliary pretreatment with ultrasound.....	18
8.4.3 Manual cleaning.....	18
8.4.4 Manual disinfection.....	18

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8.5	Assembly, inspection and care	19
8.6	Packaging systems	20
8.7	Sterilization	21
8.7.1	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	21
8.7.2	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO® V-PRO™ 1	22
8.7.3	Ethylene oxide sterilization (EO)	22
8.7.4	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®	23
8.7.5	High-Level Disinfection	23
8.7.6	CIDEX (2,4 % glutaraldehyde solution)	23
8.7.7	Resert XL HLD (2,0 % hydrogen peroxide solution)	24
8.8	Reprocessing limits	25
9	Reprocessing (for US customers only)	26
9.1	Care and handling	26
9.2	Reprocessing	27
9.2.1	Reference chart	27
9.2.2	General warnings	28
9.3	Accessories	30
9.4	Water Quality Requirements	30
9.5	Pre-cleaning	30
9.5.1	Manual cleaning	31
9.6	Assembly, inspection and care	32
9.7	Packaging systems	33
9.8	Sterilization	33
9.8.1	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	33
9.8.2	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO® V-PRO™ 1	34
9.8.3	Ethylene oxide sterilization (EO)	35
9.8.4	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E® (only valid in USA)	35
9.8.5	High-Level Disinfection	36
9.8.6	CIDEX (2,4 % glutaraldehyde solution)	36
9.8.7	Resert XL HLD (2,0 % hydrogen peroxide solution)	37
9.9	Limits of reprocessing	37

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2.2 General safety information

WARNING: Read this instruction manual thoroughly before using the equipment.

WARNING: Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.

WARNING: Perceived electric shocks may be experienced when using the camera system in conjunction with electrosurgical units.

WARNING: When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.

WARNING: Patient leakage When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. It should be noted that this is particularly important if a Type CF applied part endoscope is used, in which case a Type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

WARNING: Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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6 Operating instructions for the IMAGE1 S™ 4U camera head

6.1 Frequency of the IMAGE1 S™ 4U camera head (50/60 Hz)

The IMAGE1 S™ 4U camera head is suitable for operation at either 50 Hz or 60 Hz. The camera head must be selected to suit the local conditions.

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The IMAGE1 S™ 4U camera head has an integrated instrument coupler and integrated telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the coupler clockwise and insert the endoscope eyepiece. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the endoscope eyepiece by tightening the knurled screw on the light cable through a quarter turn.

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! **NOTE:** For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (97000213), which can be downloaded or requested at www.karlstorz.com.

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- ❗ **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.
- ❗ **CAUTION:** Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be reprocessed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.
- ❗ **CAUTION:** Chemical residue on the camera head plug can cause failure of the camera electronics.
- ❗ **CAUTION:** The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.
- ❗ **CAUTION:** The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

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- Brush: 27652

8.4 Preparation for cleaning and disinfection

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8.4.1 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened **with** disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surfaces with a dry low lint cloth.

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 **CAUTION:** V-PRO™ sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

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 **WARNING:** Medical devices must be dismantled for sterilization.

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 **CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing technique of the EO sterilizer (concentration, process control).

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**8.7.6 CIDEX (2,4 % glutaraldehyde
solution)**

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2,4 % glutaraldehyde solution for a minimum of 45 minutes at 25 °C (77 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

 **CAUTION:** Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse.

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**Reprocessing
(not validated for
US customers)**

Thorough rinsing of the camera head and the card-edge connector with sterile water is very important for removing any residual disinfectant.

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**8.7.7 Resert XL HLD (2.0 % hydrogen
peroxide solution)**

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.0 % hydrogen peroxide solution for a minimum of 8 minutes at 20 °C (68 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

! **CAUTION:** Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head with sterile water is very important for removing any residual disinfectant.

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9.2.2 General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction.

 **WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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 **CAUTION:** When preparing the cleaning and high-level disinfectant solutions, follow the manufacturer's instructions for dilution to proper concentration and temperature.

 **CAUTION:** Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices.

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- ! **CAUTION:** Inspect the camera head cable for breaks and cracks before and after preparation. Camera heads with damaged cables must NOT be processed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.
- ! **CAUTION:** Chemical residue on the camera head plug can cause failure of the camera electronics.
- ! **CAUTION:** The camera head must be thoroughly cleaned and high-level disinfected or sterilized using validated infection control procedures before use and before every subsequent application.
- ! **CAUTION:** The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.
- ! **CAUTION:** Use only one method of high-level disinfection or sterilization per reprocessing cycle.

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9.4 Water Quality Requirements

Utility water:

Water as it comes from the tap that meets specifications as defined per AAMI TIR 34:2014 Water for Reprocessing of Medical Devices. Utility water is acceptable for pre-cleaning and cleaning of the instruments.

Critical water:

Water that is extensively treated (usually by a multistep treatment process that could include a carbon bed, softening, DI and RO, or distillation) to ensure that the microorganisms and the inorganic and organic material are removed from the water; a final submicron filtration could also be part of the treatment process. Critical water is recommended for rinsing after high-level disinfection.

9.5 Pre-cleaning

Immediately after a procedure, at the point of use (e. g. operating room), wipe the camera head with a soft, lint-free disposable cloth moistened with water or diluted mild/neutral pH enzymatic cleaning solution to remove all heavy soiling.

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9.5.1 Manual cleaning

1. Thoroughly rinse the camera head using a commercial overhead sprayer nozzle with cold water to remove all gross debris. Rinse the camera head unit for a minimum of 1.5 minutes, paying special attention to the knobs, crevices, and seams. Then rinse the connecting cable for a minimum of 30 seconds.
2. Completely immerse the camera head in a diluted mild/neutral pH enzymatic cleaning solution (e.g. Enzo!). Keep immersed for at least 1 minute.
3. While immersed, use a soft bristle brush (P/N: 27652/3) to brush all the surfaces of the camera head unit. Pay close attention to the seams between the instrument coupler, focusing ring and zoom ring. Do not brush the glass window.
4. Brush for a minimum of 1 minute.
5. Then, use a low-lint wipe to wipe the entire length of the connecting cable and the glass window of the camera head.
6. Finish by brushing the entire camera head, including the entire silicon cable, card-edge connector, paying special attention to the operating buttons, screws, knobs, and grasping mechanism. Rinse the bristle of the brush if heavy debris collects on them. Keep brushing until there is no visible debris on the bristles of the brush (minimum of 1 minute).
7. After the immersion and brushing period, remove the camera head from the cleaning solution and rinse using a commercial overhead sprayer nozzle with cold water. Rinse the camera head unit for a minimum of 1.5 minutes, paying special attention to the **knobs**, crevices and seams. Then rinse the connecting cable for a minimum of 30 seconds.

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9.6 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning procedure once more.
- The device must be free from any visible signs of damage (cuts, dents, or scratches, etc.) that may trap residual debris. If the device shows any damage, contact KARL STORZ for repairs.
- Wipe the camera head connector with a cotton tip applicator soaked in 70 % isopropyl alcohol to remove any residual disinfectant. Do not allow exposed glass windows to air dry. Glass surfaces can be cleaned using a cotton tip applicator soaked in 70 % isopropyl alcohol to prevent streaks and spots.

Before sterilization and/or disinfection, coil the camera head cable into loops with a diameter of at least 6 inches (15 cm). In doing so, do not kink or twist the cable.

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9.7 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953). National regulations must be observed.

9.8 Sterilization

The following sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ 4U:

9.8.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

CAUTION: STERRAD® sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

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9.8.3 Ethylene oxide sterilization (EO)

1. Place the camera head in an FDA-cleared sterilization tray.
2. Double wrap the sterilization tray with polypropylene wrap.
3. Start sterilization under the following parameters:

Conditioning parameters

Temperature: 55 °C (131 °F)
Humidity: 70 % RH
Conditioning dwell time: 30 minutes

Sterilization parameters

Sterilant: 100 % ethylene oxide
Temperature: 55 °C (131 °F)
Gas concentration: 735 ±30 mg/l
Exposure time: 180 minutes

Aeration parameters

Time: 12 hours, minimum
Temperature: 55 °C (131 °F)

**9.8.4 Chemical low-temperature
sterilization with peracetic acid
– STERIS® system 1E®
(only valid for USA)**

1. Place the camera head in one of the following trays:
 - C1200E SS1E General Processing Tray/ Container
 - C1220E SS1E Directed Flow Processing Tray/ Container
 - C1140E SS1E Flexible Endoscope Processing Tray/Container
 - C1160E SS1E Universal Flexible Processing (Open) Tray
2. Load the tray into the SS1E and start the sterilizer.

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Reprocessing (for US customers only)

9.8.5 High-level disinfection

WARNING: High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

CAUTION: Before disinfection, the instruments must be thoroughly cleaned, rinsed, and dried.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life.

The following high-level disinfection methods have been validated and approved by KARL STORZ:

9.8.6 CIDEX (2,4 % glutaraldehyde solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.4 % glutaraldehyde solution for a minimum of 45 minutes at 25 °C (77 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

CAUTION: Do not leave the camera head in the solution for longer than 60 minutes.

To finish, immerse the camera head three times in water. Each immersion rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head and the card-edge connector with water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before inserting it into the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

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9.8.7 Resert XL HLD (2.0 % hydrogen peroxide solution)

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To finish, immerse the camera head three times in water. Each immersion rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head with water is very important for removing any residual disinfectant.

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9.9 Limits of reprocessing

The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

Always perform a functional test with the full IMAGE1 S™ system prior to using the camera head. The system should not be used if any of the accessories are damaged or non-functional.

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**Kompatible
Kamera-Kontrolleinheit**

**Compatible
camera control unit (CCU)**

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
**11 Compatible
camera control unit (CCU)**

Item	Cat. No.
IMAGE1 S™ 4U-LINK in combination with IMAGE1 S CONNECT™ or IMAGE1 S CONNECT™ II	TC 304 TC 200 TC 201

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Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Instruction Manual	96286003DER
Instruction manual IMAGE1 S™	96206286DER

* For safe and comfortable working with KARL STORZ cameras, we recommend the use of sterile drapes.

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Information on electromagnetic compatibility (EMC)

13 Information on electromagnetic compatibility (EMC)

The OCU and camera head were tested as a system.

 **WARNING:** Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Please observe the EMC instructions in this instruction manual during installation and operation.

 **WARNING:** The IMAGE1 S™ 4U camera head is only intended for use by physicians and medical assistants. This is a CISPR Class A medical equipment/system. In a domestic environment, IMAGE1 S™ 4U camera head may cause radio interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the unit or filtering the connection to the public mains network.

 **WARNING:** The use of portable and mobile RF equipment may have impact on this and other pieces of medical equipment.

 **WARNING:** The use of accessories or cables with the IMAGE1 S™ 4U camera head other than those specified in this manual may result in increased emissions or decreased immunity of the IMAGE1 S™ 4U camera head. When using accessories or cables with the IMAGE1 S™ 4U camera head other than those specified in this manual, it becomes the responsibility of the user of the IMAGE1 S™ 4U camera head to determine compliance with IEC 60601-1-2 when using this item.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic emissions		
The IMAGE1 S™ 4U camera head is intended for operation in an electromagnetic environment as defined below. The customer or user of the IMAGE1 S™ 4U camera head should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic enforcement – guidance
RF emissions acc. to CISPR 11	Group 1	The IMAGE1 S™ 4U camera head uses RF energy only for its internal function. The customer or user of the IMAGE1 S™ 4U camera head should ensure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class A	
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	The IMAGE1 S™ 4U camera head is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The IMAGE1 S™ 4U camera head is intended for use in the electromagnetic environment specified below. The user of the IMAGE1 S™ 4U camera head should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Air discharge	± 6 kV Contact discharge ± 8 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Rapid transient electrical interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	± 2 kV/1 kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV Differential mode ± 2 kV Common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ 4U camera head require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ 4U camera head be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the IMAGE1 S™ 4U camera head further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
i *NOTE: U_T is the mains alternating voltage before applying the test levels.			

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Table 4 Guidance and manufacturer's declaration – electromagnetic immunity			
The IMAGE1 S™ 4U camera head is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ 4U camera head should ensure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™ 4U camera head, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V _{rms}	$d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.			
^a The ISM (Industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. ^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges. ^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 4U camera head is used exceeds the applicable RF compliance level above, the 4U camera head should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the 4U camera head. ^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 6

Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™ 4U camera head

The IMAGE1 S™ 4U camera head is intended for use in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the IMAGE1 S™ 4U camera head can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ 4U camera head as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to < 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to < 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 818-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia Mexico S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcube@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No. 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 Mexico D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

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KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch
Room 1805-1807, Building B, 18F Beijing JFC
No. 8, Jiangmenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

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**1 Important note for users
of KARL STORZ devices
and instruments**

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1 Important note for users of KARL STORZ devices and instruments	III
2 Safety instructions	4
2.1 Explanation of warnings and cautions	4
2.2 General safety information	5
2.3 Intended Use Camera Heads	7
2.4 Contraindications	7
2.5 User qualification	7
3 Images of the equipment	8
4 Symbols employed	9
4.1 Symbols for operation	9
4.2 Label and packaging symbols	9
5 Using the D1 Camera Heads	10
5.1 Frequency of the D1 Camera Heads (50/60 Hz)	10
5.2 Using the D1 camera heads	10
5.2.1 Connecting an endoscope and light cable	10
5.2.2 Adjusting the focus	11
6 Technical description	12
6.1 Specifications	12
6.2 Storage and Operating conditions	12
6.3 Standard compliance	12
6.4 Directive compliance	12
6.5 Information on electromagnetic compatibility (EMC)	13
7 Reprocessing (not validated for US customers)	14
7.1 Care and handling	14
7.2 Reprocessing	15
7.2.1 Reference chart	15
7.2.2 General warnings	16
7.3 Accessories	18
7.4 Preparation for cleaning and disinfection	18
7.4.1 Manual wipe-down disinfection	18
7.4.2 Manual precleaning	19
7.4.2.1 Brushing the surfaces	19
7.4.2.2 Auxiliary pretreatment with ultrasound	19

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8.8.6	CIDEX (2.4 % glutaraldehyde solution)	38
8.8.7	Resert XL HLD (2.0 % hydrogen peroxide solution)	39
8.9	Reprocessing limits	39
9	Maintenance	40
9.1	Maintenance and safety check	40
9.1.1	Maintenance	40
9.1.2	Safety check	40
9.2	Servicing and repair	41
9.3	Disposal	42
9.4	Repair program	42
9.5	Limitation of liability	43
9.6	Warranty	44
10	Compatible camera control unit (CCU)	45
11	Spare parts, recommended accessories	46
12	Subsidiaries	47

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WARNING: All endoscopic devices (e.g., HF surgical probes) used in conjunction with KARL STORZ endoscopic cameras must comply with IEC 60601-2-18, the standard for medical electrical equipment – particular requirements for the safety of endoscopic equipment. The insulation of the endoscopic devices must be designed to withstand the highest HF voltage encountered during surgery. HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.

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CAUTION: When in contact with or in close proximity of an active electrosurgical instrument temporary image interference or loss can occur. If indicated by the system a hardware upgrade of TC 301 is recommended.

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2.3 Intended Use Camera Heads

The camera head is used together with the camera control unit (CCU) and is attached to a telescope during endoscopic or microscopic procedures.

2.4 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

2.5 User qualification

The D1 Camera Heads may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the D1 Camera Heads.

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
- ① "Snap-in" DCI® coupling mechanism
- ② Focusing ring/knob
- ③ User control button: for scrolling (upwards) through the menu or activating a camera function
- ④ User control button: for accessing the menu/ making a selection
- ⑤ User control button: for scrolling (downwards) through the menu or activating a camera function

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5 Using the D1 Camera Heads

5.1 Frequency of the D1 Camera Heads (50/60 Hz)

The D1 camera heads are suitable for operation at either 50 Hz or 60 Hz. The camera head must be selected to suit the local conditions.

NOTE: The camera head TH 115 is intended for operation at 50 Hz and the camera head TH 116 for operation at 60 Hz.

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The D1 camera heads feature an integrated DCI® coupling mechanism and an integrated telescope. The DCI® HOPKINS® endoscope telescopes are pushed onto the 'Snap-in' DCI® coupling mechanism. The coupler latches automatically when the telescope is inserted (see Fig. 1).

To release the telescope, press down the release mechanism ①, and pull the telescope out of the camera head (see fig. 2).

Connect special light cable 495 DC or 495 DV to the D1 camera head. To do this, push the longer end of the light cable connector into the D1 camera head until it locks into place.

Only grip the light cable by the handle piece. Never pull on the cable.

Remove the special light cable from the D1 camera head by pulling it out (see Fig. 3).

NOTE: For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (96206342ZS), which can be downloaded or requested at www.karlstorz.com. For U.S. use, please refer to PI-000011.

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7.3 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652

7.4 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

CAUTION: Inspect the camera head cable for breaks or cracks. Camera heads with damaged cables must not be sterilized or disinfected. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

7.4.1 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surfaces with a dry low lint cloth.

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Reprocessing (not validated for US customers)

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7.7.7 Resert XL HLD (2.0 % hydrogen peroxide solution)

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To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use **fresh** sterile water for each rinse. Thorough rinsing of the camera head **with** sterile water is **very** important for removing any residual disinfectant.

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**Reprocessing
(for US customers only)**

- ! **CAUTION:** Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be reprocessed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.
- ! **CAUTION:** Chemical residue on the camera head plug can cause failure of the camera electronics.
- ! **CAUTION:** The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.
- ! **CAUTION:** The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.
- ! **CAUTION:** Use only one method of high-level disinfection or sterilization per reprocessing cycle.

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8.5.1 Manual precleaning

**8.5.1.1 Auxiliary pretreatment with
ultrasound**

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8.5.2 Manual cleaning

1. Thoroughly rinse the camera head with water until all gross debris is removed.
2. Completely immerse the camera head in a diluted mild/neutral pH enzymatic cleaning solution for 1 minute.
3. While immersed, brush all surfaces of the camera head, including the card edge connector, buttons, screw, knobs and grasping mechanism with the a soft bristle brush (Art. no. 27652). Do not brush the glass window of the device, but wipe carefully using a lint-free, soft cloth. Rinse the bristles of the brush if heavy debris collects on them. Keep brushing until there is no more visible debris on the bristle of the brushes.
4. While immersed, use a lint-free, soft cloth to wipe the entire length of the cable.
5. After the immersion period, remove the camera head from the cleaning solution and rinse under running water for a minimum of 30 seconds.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged medical devices must be withdrawn from use.

! **CAUTION:** Inspect the camera head cable for breaks or cracks. Return camera heads with damaged cables to KARL STORZ for repair.

- Wipe the camera head connector with a cotton tip applicator soaked in 70 % isopropyl alcohol to remove any residual disinfectant. Do not allow exposed glass windows to air dry. Glass surfaces can be cleaned using a cotton tip applicator soaked in 70 % isopropyl alcohol to prevent streaks and spots.
- Afterwards, a functional check must be carried out.

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i **NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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8.8 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ D1 camera heads:

8.8.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

 **WARNING:** Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of dimensions and material.

 **CAUTION:** For more detailed information, please consult the user handbook of the respective device.

 **CAUTION:** STERRAD® sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

 **NOTE:** Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

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- STERRAD® 100S
 - STERRAD® NX®
‘standard’ cycle
 - STERRAD® 100NX®
‘standard’ cycle
‘DUO’ cycle
1. Place the camera head in a sterilization tray.
 2. Double wrap the tray with polypropylene wrap.
 3. Start the sterilizer using the appropriate cycle.

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CAUTION: V-PRO™ sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

1. Place the camera head in a sterilization tray.
2. Double wrap the tray with polypropylene wrap.
3. Start the sterilizer using the appropriate cycle.

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WARNING: Medical devices must be dismantled for sterilization.

NOTE: Place the camera head in a sterilization tray. Do not place any instruments on top of the camera.

CAUTION: For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing technique of the EO sterilizer (concentration, process control).

**8.8.4 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®
(only valid in USA)**

- Place the camera head in one of the following trays:
 - C1200E SS1E General Processing Tray/ Container
 - C1220E SS1E Directed Flow Processing Tray/ Container
 - C1140E SS1E Flexible Endoscope Processing Tray/ Container
 - C1160E SS1E Universal Flexible Processing (Open) Tray
- Load the tray into the SS1E and start the sterilizer.

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CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life.

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To finish, rinse the camera head three times with water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head and the card-edge connector with water is very important for removing any residual disinfectant.

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8.8.7 Resert XL HLD (2.0 % hydrogen peroxide solution)

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CAUTION: Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head with water is very important for removing any residual disinfectant.

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8.9 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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**10 Compatible camera control
unit (CCU)**

Item	Cat. No.
IMAGE1 S CONNECT™ in combination with IMAGE1 S™ X-LINK	TC200 TC301

① NOTE: The S technologies are available with IMAGE1 S™ D1 camera heads and the IMAGE1 S™.

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Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Sterile drapes* for KARL STORZ cameras	
Instruction Manual	96206564DER
Instruction manual IMAGE1 S™	96206286DER

* For safe and comfortable working with
KARL STORZ cameras, we recommend the use
of sterile drapes.

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Shawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, lower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1261, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8572 6399
E-Mail: mingwang@karlstorz.tw

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KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquirv@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F. Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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