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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SIEMENS SHENZHEN
MAGNETIC RESONANCE LTD. on the annexed DOCUMENT is
genuine.

China Council for the Promotion
of International Trade

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
证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Siemens Shenzhen Magnetic Resonance Ltd.,
Siemens MRI Center,
Gaoxin C. Ave, 2nd
Hi-Tech Industrial Park
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

User Manual Supplement for Russia

The magnetic-resonance system

MAGNETOM ESSENZA with accessories

Signature _____
Name Yu Ruo Xi
(Regulatory Affairs Specialist)

1. The name of the Medical Device.

The magnetic-resonance system MAGNETOM ESSENZA with accessories.

2. Manufacturing information of the Medical Device.

Legal manufacturer:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center,
Gaoxin C. Ave. 2nd, Hi-Tech Industrial Park,
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86 0755 26525421
www.siemens.com.cn/ssmr

Manufacturing sites:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center,
Gaoxin C. Ave. 2nd, Hi-Tech Industrial Park,
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86 0755 26525421
www.siemens.com.cn/ssmr

Siemens Healthcare GmbH
Magnetic Resonance (MR),
Henkestr. 127,
91052 Erlangen, Germany
Phone: +49 9131 84-0
www.siemens.com/healthcare

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubininskaja street 96,
115093 Moscow, Russia
Tel.: +7 495 737 1252
www.siemens-healthineers.com/ru/

3. Designation of the medical device.

3.1. *Indication:*

MAGNETOM ESSENZA system is indicated for use as a magnetic resonance diagnostic device (MRDD) that produces transverse, sagittal, coronal and oblique cross sectional images, spectroscopic images and/or spectra, and that displays the internal structure and/or function of the head, body, or extremities. Other physical parameters derived from the images and/or spectra may also be produced. Depending on the region of interest, contrast agents may be used. These images and/or spectra and the physical parameters derived from the images and/or spectra when interpreted by a trained physician yield information that may assist in diagnosis.

MAGNETOM ESSENZA system may also be used for imaging during interventional procedures when performed with MR compatible devices such as in-room displays and MR-safe biopsy needles.

3.2. *Authorized operating personnel:*

MAGNETOM ESSENZA system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training. This user training must include basics in MR technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manual.

3.3. *Contraindications and Possible side effects:*

An MR examination is contraindicated for patients with electronic or electronically conducting implants or metals, especially those containing ferromagnetic foreign bodies. These may also be patients engaged in occupations or activities which may cause accidental lodging of ferromagnetic materials or metal fragments from military activities.

Contraindications for MR examinations are:

- *Electronic implants: Pacemakers, insulin pumps*
- *Artificial heart valves; aneurysm clips*
- *Metal splinters in the eye (danger of retinal detachment)*
- *Artificial anus (anus praeter) with magnetic closure*
- *Transdermal drug patches with metallic backings*
- *Electrically conductive implants and prostheses*
- *Metallic spirals for contraception*
- *Transdermal or other similar implants for body piercings as well as magnetic piercings*

3.4. *Physiological effects:*

Due to the presence of strong alternating electromagnetic fields, patients may experience various physiological effects during MR examinations:

- Peripheral nerve stimulation through low-frequency fields of the gradient coils
- Warming of body tissue through RF fields of the RF-transmitter coil


It is generally accepted that no published evidence supporting the occurrence of cumulative and/or long-term effects after exposure to EMF emitted by the MR equipment exists.



3.5. *Conditions for use:*

The MAGNETOM ESSENZA MR system must be operated according to the indication and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training.

This user training must include basics in MR technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the Operator Manual.

3.6. *System's overview*

Note: The unit of above drawing is mm.

4. Classification of the Medical Device:

4.1. *Type of the medical device in accordance with the nomenclature for classification:*

Classification according to Annex IX, MDD

Class

Ila

Classification in accordance with Russia Nomenclature classification:

135160

4.2. *Class of potential Risk of the medical device application:*

The Medical Device is intended for transmission of energy in the form of ionizing radiation, it belong to the **Class 2a**.

4.3. *Class of the Medical Device in accordance with IEC 60601-1:*

According to IEC 60601-1:2005, the Medical Device is the **Class 1**.

4.4. *Type on contact with the human body:*

The Medical Device is non-sterile, nonexpendable.

4.5. *Software classification:*

The Medical Device's Software safety Class according to IEC 62304:2006 is **Class B**.

4.6. *Applied parts:*

According to IEC 60601-1 Cl. 7.2.10 and 8.3 the following parts are applied parts of following types:

Patient table	Type B applied part
Body coil	Type B applied part
ECG electrodes	Type BF applied part
Head Matrix Coil	Type BF applied part
Neck Matrix Coil	Type BF applied part
IsoCenter Matrix	Type BF applied part

Flex Small/Large	Type BF applied part
Head/Neck 16	Type BF applied part
Body Matrix	Type BF applied part
PA Matrix	Type BF applied part
Extremity Matrix	Type BF applied part
14-ch Extremity Coil	Type BF applied part
8 Ch Foot Ankle	Type BF applied part
8 Ch Wrist	Type BF applied part
Shoulder Array Small/Shoulder Array Large	Type BF applied part
Breast Matrix	Type BF applied part
Endorectal	Type BF applied part
Flex Large	Type BF applied part
Flex Small	Type BF applied part
4 Ch Special Purpose	Type BF applied part
Physiological measurement unit:	
- ECG sensor	Type BF applied part
- Respiratory sensor	Type BF applied part
- Pulse sensor	Type BF applied part

The system is not classified as defibrillator-proof applied part (IEC 60601-1:2012 Clause 8.5.5).

4.7. IP protection class:

IP protection class according to IEC 60529 >>> IP X0

4.8. Laser positioner:

Laser light Class: **2M** (according to IEC 60825 - 1:2007).

4.9. Operating mode:

Operating conditions according to IEC 60601-1.

System: Continuous operation mode

Patient table: 20 cm/s (Max. speed)

5. Technical Data: Specifications and Characteristics.

Tolerances: $\pm 5\%$ (if not indicated another information)

5.1. Technical Data:

Magnet system

- Short bore, patient-friendly design with flared opening
- Simple design and active shielding for easy installation
- Integrated external interference shielding provides permanent protection against potential interferences through ferromagnetic objects moved into the proximity of the magnet
- Operating mode: Standard operating mode according to IEC 60601-2-33

Magnet bore	
Magnet length	131 cm
Inner diameter*	60 cm
Inner tunnel length	103 cm
Length of flared opening	28 cm
Magnet bore diameter**	82 cm
System length	147 cm

* incl. shim coils, gradient coil, RF body coil

**without shim coils, gradient coil, RF body coil

Magnet parameters	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	< 0.1 ppm/h

Homogeneity		
Volume	Guaranteed	Typical
10 cm DSV	≤ 0.2 ppm	0.04 ppm
20 cm DSV	≤ 0.22 ppm	0.16 ppm

30 cm DSV	≤ 1.2 ppm	1.1 ppm
40 × 40 × 30 cm ³	≤ 2.2 ppm	1.9 ppm
45 × 45 × 30 cm ³	≤ 7.5 ppm	6.0 ppm
Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).		
40 × 40 × 40 cm ³	≤ 5 ppm	

DSV = Diameter spherical volume

3D Shim	
■ Patient-specific automated shim ■ For optimization of the magnetic field homogeneity	
Time to shim	Approx. 22s

Active Shielding (AS)		
Fringe field	0.5 mT*	4.0 m × 2.5 m
(axial × radial)	0.1 mT	5.8 m × 3.4 m

* pacemaker safety limit

External Interference Shield (E.I.S.)
Patented shielding system integrated into the magnet
Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects)

Magnet cooling system	
Refill interval (typical)	10 years*
Boil-off rate (typical)	0.0 l/year*
Helium capacity	approx. 1400 liters
Minimum helium level	40%
Stainless steel	Cryostate

* for typical clinical use, depending on sequences and operating time

5.2. Patient handling:

The patient table is available in two different configurations:

- Matrix configuration (with vertical drive)
- Matrix configuration (without vertical drive)

The respective configuration is set by Siemens Service during MR system installation.

Patient table	
Max. patient weight (including vertical movement)	200 kg
Max. scan range	140 cm
Vertical table movement*	
Range	55 × 89 cm
Max. vertical movement time	13 s
Longitudinal table movement	
Max. range	205.5 cm
Max. speed	20 cm/s
Positioning accuracy	± 0.8 mm

* optional

Cushions and positioning aids

Table cushion 1	
Manufacturer	Leidel & Kracht
PN	10168111
Weight	1.15kg
Dimensions	445.5 mm × 494 mm × 20 mm
Table cushion 2	
Manufacturer	Leidel & Kracht
PN	10168113
Weight	2.0kg
Dimensions	1180 mm × 494 mm × 20 mm
Table cushion 3	
Manufacturer	Leidel & Kracht
PN	10168112

Set of

Weight	2.0kg
Dimensions	523 mm × 494 mm × 20 mm
Table cushion 4	
Manufacturer	SIEMENS
PN	07580231
Weight	1.8kg
Dimensions	1100 mm × 490 mm × 20 mm
2 belt straps;	
Manufacturer	Siemens Medical Equipment Ltd. Shanghai
PN	10167737
Weight	0.3kg
Dimensions	1000mm×50 mm ×1 mm
2 different-sized sandbags, 2 each	
Manufacturer	SIEMENS
PN	05516104
Weight	2.5kg; 1.5kg
Dimensions	150 mm ×400mm; 150 mm ×300mm
Intermediate connector	
Manufacturer	SIEMENS
PN	07573939
Weight	0.1kg
Dimensions	50mm×40 mm ×15 mm

positioning aids

1 knee wedge	
Manufacturer	PolyForm GmbH
PN	10185620
Weight	1.25kg
Dimensions	557 mm × 395 mm × 165 mm
1 positioning roll	
Manufacturer	PolyForm GmbH
PN	10835532
Weight	1.4kg
Dimensions	220 mm × 220 mm × 399 mm
1 back-of-the head support	
Manufacturer	PolyForm GmbH
PN	10185779

Weight	0.45kg
Dimensions	250 mm × 250 mm × 10 mm
2 large triangles;	
Manufacturer	PolyForm GmbH
PN	10835530
Weight	0.4kg
Dimensions	235 mm × 176 mm × 117 mm
2 small triangles	
Manufacturer	PolyForm GmbH
PN	10606835
Weight	0.4kg
Dimensions	176 mm × 94 mm × 48 mm
1 wedge;	
Manufacturer	PolyForm GmbH
PN	10835531
Weight	0.35kg
Dimensions	237 mm × 176 mm × 67 mm
1 lordosis cushion	
Manufacturer	PolyForm GmbH
PN	10185780
Weight	0.4kg
Dimensions	345 mm × 247 mm × 40 mm

Head Matrix accessories Coil

Relaxation cushion	
Manufacturer	PolyForm GmbH
PN	10185618
Weight	0.4kg
Dimensions	310 mm × 250 mm × 80 mm
Head stabilization wedge	
Manufacturer	PolyForm GmbH
PN	10185619
Weight	0.05KG
Dimensions	140 mm × 80 mm × 75 mm
Head support.	
Manufacturer	PolyForm GmbH
PN	08624731
Weight	0.1KG

Dimensions	200 mm ×60 mm ×200 mm
------------	-----------------------

Neck Matrix Coil accessories

Neck cushion	
Manufacturer	PolyForm GmbH
PN	08624723
Weight	0.1KG
Dimensions	290 mm ×213 mm ×60 mm
2 positioning aids for the neck	
Manufacturer	PolyForm GmbH
PN	08625472
Weight	0.1KG
Dimensions	155 mm ×50 mm ×50mm

16-channel Head/Neck Coil accessories

Relaxation cushion	
Manufacturer	PolyForm GmbH
PN	10606510
Weight	0.3KG
Dimensions	230 mm ×275 mm ×155 mm
Head/Neck cushion	
Manufacturer	PolyForm GmbH
PN	08624723
Weight	0.1KG
Dimensions	290 mm ×213 mm ×60 mm
Head stabilization wedges	
Manufacturer	PolyForm GmbH
PN	07101822
Weight	0.03KG
Dimensions	135 mm ×60 mm ×50 mm
Clamp cushions	
Manufacturer	PolyForm GmbH
PN	08626314
Weight	0.2KG

Dimensions	120 mm ×60 mm ×70 mm
Mirror carrier with mirror Head/Neck 16	
Manufacturer	SIEMENS
PN	10352472
Weight	0.25KG
Dimensions	186 mm ×110 mm ×120 mm

Body Matrix Coil accessories

Mirror carrier with mirror Head/Neck 16	
Manufacturer	Pi-Products
PN	07580900
Weight	0.4KG
Dimensions	509 mm ×311 mm ×14 mm
6 belts	
Manufacturer	Siemens Shanghai Medical Equipment Ltd.
PN	10168278
Weight	0.25KG
Dimensions	1000mm×50 mm ×1 mm
2 straps.	
Manufacturer	Siemens Shanghai Medical Equipment Ltd.
PN	05514232
Weight	0.11KG
Dimensions	2000mm×50 mm ×1 mm

Breast Matrix Coil accessories

2 large compression attachments	
Manufacturer	PolyForm GmbH.
PN	05513192
Weight	0.25KG
Dimensions	138mm×95 mm ×36 mm
2 small compression attachments	
Manufacturer	PolyForm GmbH.
PN	05514448
Weight	0.06KG

Dimensions	122mm×50 mm ×7 mm
Positioning aids;	
Manufacturer	PolyForm GmbH.
PN	05513192
Weight	0.2KG
Dimensions	138mm×95 mm ×30 mm
Soft pad	
Manufacturer	SIEMENS
PN	05513184
Weight	0.55KG
Dimensions	500mm×490 mm ×10 mm
Head cushion	
Manufacturer	Pi Products
PN	08625480
Weight	1.31KG
Dimensions	330 mm × 302 mm × 85 mm
Wedge cushion (abdominal support);	
Manufacturer	PolyForm GmbH
PN	10165048
Weight	2.45KG
Dimensions	487 mm × 347 mm × 157 mm
Special head support	
Manufacturer	PolyForm GmbH
PN	10351843
Weight	0.85KG
Dimensions	420 mm × 200 mm × 40 mm

PA Matrix Coil accessories

Belt pad	
Manufacturer	PolyForm GmbH
PN	10164888
Weight	4KG
Dimensions	835 mm ×295 mm ×180 mm
2 knee cushions	
Manufacturer	PolyForm GmbH
PN	07098168
Weight	0.26KG

Dimensions	155 mm × 115 mm × 77 mm
4 flat wedges;	
Manufacturer	PolyForm GmbH
PN	07098150
Weight	0.4KG
Dimensions	250 mm × 150 mm × 51 mm
2 flat cushions.	
Manufacturer	PolyForm GmbH
PN	07098176
Weight	0.4KG
Dimensions	300 mm × 100 mm × 20mm

Focus Shoulder Array Coil accessories

Cushion shoulder.	
Manufacturer	PolyForm GmbH
PN	10351842
Weight	0.08KG
Dimensions	325 mm × 284 mm × 15 mm
Special head support	
Manufacturer	PolyForm GmbH
PN	10351843
Weight	0.3KG
Dimensions	420 mm × 200 mm × 40 mm
Shoulder support cushion	
Manufacturer	PolyForm GmbH
PN	10352369
Weight	0.05KG
Dimensions	364 mm × 220 mm × 10 mm

Extremity Matrix Coil accessories

Coil pad, 0.25	
Manufacturer	PolyForm GmbH
PN	10351926
Weight	0.15KG
Dimensions	250 mm × 200 mm × 5 mm
Coil pad, 0.50	
Manufacturer	PolyForm GmbH

PN	10351927
Weight	0.25KG
Dimensions	250 mm × 200 mm × 10 mm
Ankle pad	
Manufacturer	PolyForm GmbH
PN	10351923
Weight	0.45KG
Dimensions	200 mm × 150 mm × 220 mm
Ramp Pad	
Manufacturer	PolyForm GmbH
PN	10351928
Weight	0.65KG
Dimensions	200 mm × 187 mm × 260 mm

14-channel Extremity Coil accessories

Coil pad, 0.5	
Manufacturer	PolyForm GmbH
PN	10351926
Weight	0.15KG
Dimensions	250 mm × 200 mm × 5 mm
Coil pad, 1.0	
Manufacturer	PolyForm GmbH
PN	10351927
Weight	0.25KG
Dimensions	250 mm × 200 mm × 10 mm
Knee cushion	
Manufacturer	PolyForm GmbH
PN	10500140
Weight	0.6KG
Dimensions	180 mm × 160 mm × 220 mm
Air cushion	
Manufacturer	SIEMENS
PN	01597157
Weight	0.3KG
Dimensions	300 mm × 190 mm × 50 mm
Ramp pad.	
Manufacturer	PolyForm GmbH
PN	10353257

Weight	0.15KG
Dimensions	205mm×200mm×137mm

4-channel Flex Coil, large accessories

4-channel Flex Coil, small accessories

Hand Wedge;	
Manufacturer	PolyForm GmbH
PN	IPAT cushion flex large:10185345 IPAT cushion flex small:10185346
Weight	IPAT cushion flex large: 1.0kg IPAT cushion flex small: 0.3kg
Dimensions	IPAT cushion flex large: 504mm×212mm×14mm IPAT cushion flex small: 354mm×162mm×14mm
2 fleece straps	
Manufacturer	SIEMENS
PN	05514232
Weight	0.11KG
Dimensions	2000mm×50 mm ×1 mm

Endorectal Coil Interface accessories

Endo cushion	
Manufacturer	Bayer Medical Care, Inc
PN	07582393
Weight	1.9KG
Dimensions	450 mm × 160 mm × 130 mm

8-channel Foot/ Ankle Coil accessories

Knee cushion	
Manufacturer	PolyForm GmbH
PN	10353257
Weight	0.15KG
Dimensions	205mm×200mm×137mm
Heel cushion	

Manufacturer	PolyForm GmbH
PN	10353290
Weight	0.05KG
Dimensions	155mm×102mm×125mm
Air cushion	
Manufacturer	SIEMENS
PN	01597157
Weight	0.15KG
Dimensions	300mm×190mm×50mm

8-channel Wrist Coil accessories

Wrist pad	
Manufacturer	PolyForm GmbH
PN	10352471
Weight	0.05KG
Dimensions	125 mm ×46 mm ×85 mm
Palm pad	
Manufacturer	PolyForm GmbH
PN	10352471
Weight	0.02KG
Dimensions	110 mm ×160 mm ×7 mm
Upper pad	
Manufacturer	PolyForm GmbH
PN	10352520
Weight	0.02KG
Dimensions	110 mm ×160 mm ×7 mm
Air cushion	
Manufacturer	SIEMENS
PN	01597157
Weight	0.15KG
Dimensions	300 mm ×190 mm ×56 mm

Coil Storage Cart

Manufacturer	SIEMENS
PN	10606590
Weight	150KG
Dimensions	1520 mm × 730 mm × 1400 mm
Wheels	4 mm

Patient Supervision TV

Manufacturer	10164365: Shenzhen Skyworth Qunxin Security Technology Co., Ltd 10433479: CONRAC GmbH
Model	10164365: M17LP 10433479: 6012-MD003-001
Weight	10164365: 5.3kg 10433479: 4.0kg
Dimensions	10164365: 384.2 mm x 420.4 mm x 180mm 10433479: 277.7 mm x 294.3 mm x 142.6mm

Patient comfort	
Inner diameter of magnet bore*	60 cm
In-bore lighting	can be set at 3 different levels
In-bore ventilation	can be set at 3 different levels
In-bore intercom	including loudspeaker, microphone and earphones

* incl. shim coils, gradient coil, RF body coil

5.3. Hearing protection data:

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4 -2010 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) mode.

Patient noise

Patients require hearing protection with SNR = 19 dB or more.

Noise for personal in the examination room 0

Noise measured: 87.1 dB(A).

Hearing protection: SNR = 6 dB or more.

5.4. Gradient system:

General features	
Actively shielded (AS), super-compact whole-body gradient coil system	
Extremely low eddy currents	
Shortest TE and TR for high image contrast and best signal-to-noise ratio	
Water-cooled coil and amplifier for maximum performance	
Gradient coil length	123.5 cm

Patient-specific shimming
Active shim via 3 linear channels
Includes non-linear 2nd order channels for advanced order shim

Gradient amplifier		
- Water-cooled, highly compact, modular design - Ultra-fast solid-state technology with very low switching losses		
V-engine	Max. output voltage*	1350 V
	Max. output current*	150 A

* Values for each of the 3 gradient axes

Gradients (Tim 25 × 8)(Tim 46 × 16) V-engine	
Performance	
Max. amplitude	30 mT/m
Min. rise time	300 µs
Max. slew rate	100 T/m/s
Vector gradient performance (vector summation of all 3 gradient axes)	
Max. eff. amplitude	52 mT/m
Max. eff. slew rate	173 T/m/s
Gradient duty cycle	100%

5.5. RF system:

Frequency stability (5 min)	$\pm 2 \times 10^{-10}$
Frequency control	32 bits (0.015 Hz)
Phase control	16 bits (0.006 degrees)

Transmitter path	
Transmit amplitude	16 bit control 50 ns resolution
Gain stability (5 min)	<0.15 dB (low signal path) <0.40 dB (power amplifier)

Transmit amplifier	
Extremely compact water-cooled	
Bandwidth	500 kHz
Peak power	15 kW

Receive path	
Bandwidth	500 Hz– 1 MHz (for each channel)
Receiver signal resolution	32 bits
Number of independent receiver channels	8/16
Sampling rate (Sampling resolution)	10 MHz (100 ns)
Preamplifier noise figure	<0.55 dB
Total gain	35 dB/50 dB with automatic control
Dynamic range	162 dB instantaneous at receiver 169 dB with automatic control

5.6. RF coils:

Matrix coils (standard)	specifications
Head Matrix	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10131483
RF field	63.6MHz
Dimensions	359 mm x 385 mm x 307 mm(L x W x H)
Weight	4.6 kg
Applications	● Head examinations ● MR Angiography ● Combined head/neck examinations ● TMJ (temporo mandibular joints)
16-Channel Head/Neck coil	
Applications	● Head examinations ● Cervical spine ● Neck ● Larynx/Esophagus ● MR Angiography ● Mediastinum ● Combined head/neck examinations ● TMJ (temporo mandibular joints)
Neck Matrix	
Applications	● Cervical spine ● Neck ● Larynx/Esophagus ● MR Angiography ● Mediastinum ● Combined head/neck examinations
IsoCenter Matrix	

Applications	Imaging of multiple body parts in single and multi-step examinations in matrix coil image quality
--------------	---

Flexible coils (standard)	specifications
The Flex coil interface.	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10498864
Weight	2.05kg
Dimensions	865 mm × 83 mm × 43mm (L × W × H)
4-Channel Flex Coil, large	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd.
Model	10498865
RF field	63.6MHz
Weight	550 g
Dimensions	516 mm x 224 mm (L × W)
Applications	Imaging of large regions such as hip and knee
4-Channel Flex Coil, small	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd.
Model	10498866
RF field	63.6MHz
Weight	450 g
Dimensions	366 mm x 174 mm(L × W)
Applications	Imaging of small regions such as wrist, elbow and ankle

Matrix coils (optional)	specifications
Body Matrix	

Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10131487
RF field	63.6MHz
Weight	1.5 kg; "patient-felt" weight of coil only 950 g
Dimensions	322 mm × 520 mm × 40 mm (L × W × H)
Applications	● Thorax (incl. heart) ● Abdomen ● Pelvis ● Hip
PA Matrix (Peripheral Angio Matrix)	
Manufacturer	Siemens Healthcare GmbH
Model	07579910
RF field	63.6MHz
Weight	5.75 kg
Dimensions	970 mm × 300–600 mm × 270 mm (L × W × H)
Applications	● High resolution angiography of both legs incl. pelvis with highest signal-to-noise ratio ● Visualization of the iliac arteries and aorta ● Bilateral examinations of long bones of the legs
Breast Matrix	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10131496
RF field	63.6MHz
Weight	7.5 kg
Dimensions	500 mm × 520 mm × 135 mm (L × W × H)

Applications	• Simultaneous high resolution imaging of both breasts in all directions • Uni- or bi-lateral imaging of the breasts in sagittal direction • Visualization of axilla region
Extremity Matrix	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10351919
RF field (MHz)	63.6MHz
Weight	6.5 kg
Dimensions (mm)	536 mm × 240 mm × 297 mm (L × W × H)
Applications	• Knee • Ankle • Peripheral MR Angiography

Dedicated coils (optional)	specifications
Focus Shoulder Array, large/small	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	Focus Shoulder Array coil, large 10351766 Focus Shoulder Array coil, small 10352367
RF field	63.6MHz
Weight	Focus Shoulder Array coil, large 2.2 kg (Base Plate Weight: 2.8 kg) Focus Shoulder Array coil, small 2.1 kg (Base Plate Weight: 2.8 kg)
Dimensions	525 mm × 200 mm (L × W)

Dedicated coils (optional)	specifications
Applications	● Best visualization of small anatomical structures (e.g. labrum) ● Higher SNR and better field homogeneity ● Reduced slice thickness and measurement times
Endorectal Coil	
Manufacturer	Endo Interface: Bayer Medical Care Inc. Endo adapter: Siemens Shenzhen Magnetic Resonance Ltd.
PN	Endo Interface: 07582385 Endo adapter: 10151822
RF field (MHz)	63.6MHz
Weight	Endo Interface: 0.4 kg Endo adapter: 0.7kg
Dimensions (mm)	Endo Interface: 100mm × 67 mm × 34mm Endo adapter: 123 mm × 63 mm × 34mm
Applications	● Excellent visualization of the prostate, colon, rectum and cervix without the risk of invasive procedures ● Non-invasive preoperative diagnostic evaluation and treatment planning
Storage and transportation requirements:	
Temperature	-10°C to +50°C
Pressure	70 kPa to 106 kPa
Humidity	5% to 95% RH
8-Channel Foot-Ankle Coil	

Dedicated coils (optional)	specifications
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10498471
RF field (MHz)	63.6MHz
Weight	3.2 kg
Dimensions (mm)	315 mm × 213 mm × 286 mm (L × W × H)
Applications	● High resolution foot and ankle imaging
8-Channel Wrist coil	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd.
Model	10352487
RF field	63.6MHz
Weight	2.1 kg
Dimensions	280 mm × 189 mm × 126 mm (L × W × H)
Applications	● High resolution wrist and hand imaging
4-Channel Special Purpose Coil	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd.
Model	10499077
RF field	63.6MHz
Weight	300 g
Dimensions	418 mm × 127 mm × 52 mm (L × W × H)
Applications	● Excellent visualization of carotid plaque ingredients with very high signal to noise ratio
2/4-Channel Sentinelle Breast Coil	
Manufacturer	Sentinelle Medical Inc
Model	10499208
RF field	63.6MHz

Dedicated coils (optional)	specifications
Weight	10 kg
Dimensions	870 mm × 550 mm × 145 mm (L × W × H)
Applications	● High-resolution 2D and 3D imaging of both breasts in all directions ● Breast biopsy with medial and lateral
Grid of Sentinelle BreastCoil	
Manufacturer	Invivo Corporation
PN	10499215
Shelf life:	single use only
Weight	0.5kg
Dimensions	193.421 mm × 127 mm × 11.1252 mm (L × W × H)
Storage and transportation requirements	
Temperature	-10°C to +50°C
Pressure	70 kPa to 106 kPa
Humidity	5% to 95% RH
14-Channel Extremity Coil	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd.
Model	10353300
RF field	63.6MHz
Weight	7.5 kg
Dimensions	536 mm x 240 mm x 313 mm (L × W × H)
Applications	● High SNR of the coil supports high resolution imaging of knee, foot, ankle, hand and wrist
8-channel Foot/ Ankle Coil	
Manufacturer	Siemens Shenzhen Magnetic Resonance Ltd
Model	10498471

Dedicated coils (optional)	specifications
RF field	63.6MHz
Weight	3.2 kg
Dimensions	315 mm × 213 mm × 286 mm (L × W × H)
Applications	High resolution foot and ankle imaging

5.7. Sequences:

Spin Echo family of sequences

- Spin Echo (SE)—Single, Double and Multi Echo (up to 32 echoes); Inversion Recovery (IR)
- 2D/3D Turbo Spin Echo (TSE)—Restore technique for shorter TR times while maintaining excellent T2 contrast; TurboIR: Inversion Recovery for STIR, DarkFluid T1 and T2, TrueIR; Echo Sharing for dual-contrast TSE
- 2D/3D HASTE (Half-Fourier Acquisition with Single Shot Turbo Spin Echo)— Inversion Recovery for STIR and DarkFluid contrast
- SPACE (Sampling Perfection with Application Optimized Contrasts using different flip angle evolutions) — isotropic 3D T2 or PD weighted imaging acquired with PD var, T2 var or constant mode.
- TSE-based 2-point DIXON — fat and water separation in the body regions where desirable spectral fat suppression is difficult to achieve.

Gradient Echo family of sequences 0

- 2D/3D FLASH (spoiled GRE)—dual echo for in-/opposed phase imaging
- 3D VIBE (Volume Interpolated Breath-hold Examination)—quick fat saturation; double echo for in-phase/opposed phase 3D imaging; Inline Breast Evaluation
- 2D/3D MEDIC (Multi Echo Data Image Combination) for high resolution T2 weighted orthopedic imaging and excellent contrast
- 2D/3D TurboFLASH—3D MPRAGE; single shot T1 weighted imaging e.g. for abdominal imaging during free breathing
- 3D GRE for field mapping
- 2D/3D FISP (Fast Imaging with Steady State Precession)
- 2D/3D PSIF—PSIF Diffusion
- Echo Planar Imaging (EPI)—diffusion weighted; single shot SE and FID e.g. for BOLD imaging and Perfusion weighted imaging;

- 2D/3D TrueFISP—Shared Phases Real-time TrueFISP (w/o ECG gating)
- ce-MRA sequence with Inline subtraction and Inline MIP
- 2D/3D Time-of-Flight (ToF) Angiography—single slab and multi slab; triggered and segmented*
- 2D/3D Phase Contrast Angiography—multiple velocity encoding; triggered
- CV — truFISP, FLASH and turboFLASH could be derived from it for cardiac imaging.*

* for triggering the PMU wireless Physio Control option is required

5.8. Computer system:

Host computer (syngo Acquisition Workplace)	
New state-of-the-art CPU generation	
Intuitive Windows interface supports multi-tasking for fast operation	
Processor	Quad Core Xeon™
Manufacturer	HP
Model	Z440
Clock rate	≥ 3.6 GHz
Main memory (RAM)	≥ 32 GB
1 st hard disk (system SW)	≥ 300 GB
2 nd hard disk (data base)	≥ 300 GB
3 rd hard disk (images)	300 GB (approx. 440000 images, matrix 256 ² and 512 ² uncompressed)
CD-R writer	approx. 4000 images 256 ² ; DICOM Standard, ISO 9660
DVD-R writer	approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
Media drives	CD/DVD-ROM drive
Color LCD monitor	Screen matrix 1280 × 1024 pixels
Manufacturer	Shenyang Torch-Bigtide Digital Technology Co., Ltd.
Model	19" LCD Monitor: DSC1913-D

Measurement and Reconstruction System (MARS)
System configuration with 8 receiver channels

CPU	1 x Quad Core Xeon™
Clock rate	≥ 3.2 GHz
Main memory (RAM)	≥ 16 GB
1 st hard disk (system SW)	≥ 500 GB
2 nd hard disk (data base)	≥ 500 GB
Reconstruction speed	11800 recons per second (256 ² PeFT, full FoV) 54140 recons per second (256 ² PeFT, 25% FoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
System configuration with 16 receiver channels	
CPU	Quad Core Xeon™
Clock rate	≥ 3.2 GHz
Main memory (RAM)	≥ 16 GB
1 st hard disk (system SW)	≥ 500 GB
2 nd hard disk (data base)	≥ 500 GB*2
Reconstruction speed	11800 recons per second (256 ² PeFT, full FoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets

syngo MR Workplace (optional)	
New state-of-the-art CPU generation	
Intuitive Windows interface supports multi-tasking for fast operation	
Processor	Intel Xeon™ W3520
Manufacturer	HP
Model	Z440
Clock rate	≥ 3.6 GHz
Main memory (RAM)	≥ 8 GB
1 st hard disk (system SW)	≥ 300 GB
2 nd hard disk (images)	≥ 300 GB
CD-R writer	approx. 4000 images 256 ² ; DICOM Standard, ISO 9660

syngo MR Workplace (optional)	
DVD-R writer	approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
Media drives	CD/DVD-ROM drive
Color LCD monitor	Screen matrix 1280 × 1024 pixels
Manufacturer	Shenyang Torch-Bigtide Digital Technology Co., Ltd.
Model	19" LCD Monitor: DSC1913-D

5.9. Network/data coupling:

The MR system provides the level of safety according IEC 60950-1 outside patient environment. All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

5.10. Image viewing and filming:

Image recall

- Images are stored in DICOM format, allowing fast image access and recall
- Combines images in a new series, e.g. post-processed images
- Loads images automatically
- Inline Display

Image display

- Various display layouts selectable
- Up to 3 patients can be simultaneously active in the viewer
- Image annotation and labeling
- Non-interpolated display
- Fast paging through up to 500 images with 15 images/s for full screen display

Windowing

- Freely selectable window width and center
- Windowing on succeeding images
- Auto-windowing for optimized contrast

- Saves and sends window values

Automatic movie for cine display

Interactive movie paging by dragging the mouse or automatic movie mode by clicking the icon

Evaluation

Parallel evaluation of up to 40 regions of interest

- Circle
- Rectangle
- Freehand ROI
- Pixel lens with position marker
- Statistical evaluation
- Area
- Standard deviation
- Mean value
- Min/max values
- Image scrolling
- Magnification
- Distance
- Angle

2D Post-processing

Image manipulations

- Reversal of gray-scale values
- Image rotation by 90° or by user-defined angle
- Flip horizontally/vertically
- Image zoom and pan
- Shutter
- Annotation

Position display

Displays measured slice positions on localizer image and selected series

Argus viewer

Viewing software for cardiac MR studies and large data sets

- Efficient cine review of cardiac and other dynamic data sets
- Multiple sorting options
- Single movie as well as 2, 4, or 8 simultaneous slices together in movie mode
- Rapid avi creation of 1 to 8 slices simultaneously
- Creates and edits DICOM structured reports

Dynamic analysis

Arithmetic operations on images and series (e.g. for evaluation of contrast media studies)

- Addition, subtraction, multiplication, division of single images and whole series
- Arithmetic mean and standard deviation across a range of selected images
- Calculation of T1 and T2, and logarithmic images
- Differentiation/integration of selected images
- Calculation of a mean slope image from a range of selected images
- Calculation of z-score (t-test) images for evaluation of BOLD imaging data (Blood oxygenation level Dependent)
- Time-to-peak evaluation (TTP)
- ADC maps

Several evaluation functions may be started consecutively in the background

Image Filter

The Image Filter software uses adaptive filtering and adjusts to the local image content automatically

Mean Curve

Time-intensity analysis for contrast-enhanced examinations

- Creates and edits DICOM structured reports

Filming

- Connection via DICOM Basic Print
- Interactive filming
- Filming parallel to other activities
- Independent scanning and documentation—no wait time due to camera delays
- Freely selectable positioning of images onto virtual film sheet

- Selectable various film layouts
- Mother-in-Child display
- Windowing, image zoom and pan on film sheet
- Configurable image text
- Simultaneous handling of multiple film jobs
- Up to 100 virtual film sheets

Printing on paper

Interface and software for printing images on paper (laser printer not included)

Grey levels and color printing supported

Data format Postscript Level 2

5.11. 3D Post-processing:

MPR—Multiplanar Reconstruction 0

Real-time multiplanar reformatting of secondary views

- Viewing perspectives: sagittal; coronal; axial; oblique; double oblique; curved (freehand)
- Reconstruction along polygon and/or curved (freehand) cut lines
- Reconstruction based on reconstructed planes possible
- Reconstruction of user-defined ranges of parallel, radial or freehand cuts
- Selectable slice thickness and slice increment of reconstructed images
- Storing of post-processing protocols
- Annotations and 2D evaluations such as distance and ROI

MIP—Maximum Intensity Projection

3D reconstructions of vessels from a 3D data set, or a 2D sequential slice data set (acquired with dedicated MR Angiography sequences)

- Volume of Interest (Vol) defined to increase reconstruction speed and to improve image quality
- Freehand MIP: virtually eliminating interfering overlays by manually drawn, irregular contours in all dimensions
- Arbitrary views along any direction defined interactively with mouse-driven virtual trackball
- Multiple view angles around any orthogonal axis
- Projections displayed as single images, as interactive movie or by fast paging
- MIP thin/MIP thick

MinIP—Minimum Intensity Projection

Similar to MIP but reconstructs the minimum intensity (e.g. for Dark Blood techniques)

SSD—Shaded Surface Display

Three-dimensional display of surfaces, such as contrast-enhanced vessels

- Selectable variable threshold values
- Multiple view angles around any orthogonal axis
- Rectangular and irregular Volumes of Interest (Vol) can be defined to improve image quality

5.12. Ambient conditions for use:

Control room	
Temperature	15–30° C
Rel. humidity	40 to 80%
Water content	11.5g/kg
Pressure	700 – 1060 hPa

Electronics room	
Temperature	15–30° C (with integrated cabinet temperature control)
Rel. humidity	40 to 80%
Water content	11.5g/kg
Pressure	700 – 1060 hPa

Examination room	
Temperature	18–24° C
Rel. humidity	40 to 60%
Water content	11.0g/kg
Pressure	700 – 1060 hPa

5.13. Specific absorption rates:

The specific absorption rate for the RF power and the rate of change of the gradient fields are checked according to the requirements of norm IEC 60601-2-33.

5.14. Cooling system:

Water consumption	> 30 l/min*
Heat dissipation to water	≤ 21 kW

* Water temperature: 6–12° C (43–54 °F)

5.15. Line power supply:

Variable	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	± 10%
Frequency	50 Hz, 60 Hz	± 1 Hz
Gross power	50 kVA	
On-site protection	at least 140 A	
Short-time	50 kVA	
Long-time	22 kVA	

For all products, line power has to be supplied via an on-site system contactor or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710.¹

1 In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein—to the extent permitted by relevant local and national laws—in order to ensure the safety of operating personnel, patients, and third parties.

Gross power for China only: 45KVA is only for 380 V/50 Hz

Power rating plate

380V/400V/420V/440V/460V/480V~	
3PHASE	50/60 Hz
SHORT-TIME	kVA:50
LONG-TIME	kVA:22

5.16. Power and energy consumption:

System off	≤ 6.5 kW
Stand by	≤ 6.7 kW
Ready for measurement	≤ 9.0 kW

Typical examination	≤ 9.9 kW
---------------------	---------------

5.17. General classifications:

Protection class	I
Type of applied parts	
Patient table	Type B applied part
Body coil	Type B applied part
ECG electrodes	Type BF applied part
Head Matrix Coil	Type BF applied part
Neck Matrix Coil	Type BF applied part
IsoCenter Matrix	Type BF applied part
Flex Small/Large	Type BF applied part
Head/Neck 16	Type BF applied part
Body Matrix	Type BF applied part
PA Matrix	Type BF applied part
Extremity Matrix	Type BF applied part
14-ch Extremity Coil	Type BF applied part
8 Ch Foot Ankle	Type BF applied part
8 Ch Wrist	Type BF applied part
Shoulder Array Small/Shoulder Array Large	Type BF applied part
Breast Matrix	Type BF applied part
Endorectal	Type BF applied part
Flex Large	Type BF applied part
Flex Small	Type BF applied part
4 Ch Special Purpose	Type BF applied part
Physiological measurement unit:	
- ECG sensor	Type BF applied part
- Respiratory sensor	Type BF applied part
- Pulse sensor	Type BF applied part

IP protection class according to IEC 60529	IP X0
Explosion protection	The MR system is not intended for operation in areas prone to explosion (e.g., highly flammable mixtures of anaesthesia gases with air or oxygen or nitrous oxide)
Operating mode	Continuous operation mode
Disinfection Receptors/ system components	Disinfectants without alcohol, ether
Degree of safety in the presence of a flammable anesthetics mixture with air or with oxygen or with nitrous oxide	No AP or APG category equipment
Mains operated equipment with additional power sources	None

5.18. Dimensions:

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipation [kW]
Magnet 1.5 Tesla AS (incl. Helium)	211.4	131	214	3150	
Magnet in operation, incl. gradient coil, body coil, patient table, and covers	231	411	218	4365	3
Patient table	55.5 (top)	210.5	55 - 89	270	
Required min. room height clearance			240		
Min. transport dimensions	226*	183*	216*		
syngo Acquisition Workplace (table + monitor)	120	80	117 (72 + 45)		2
Host computer	17	45	46	21.5	
syngo MR Workplace (table + monitor)	120	80	117 (72 + 45)		
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	95	65	189**	700	3
Active cooling system	100	65	189**	435***	
Passive cooling system	80	65	189**	280***	

* without electronics, table and table bracket

** without attachments

*** without helium compressor

5.19. *Software:*

MAGNETOM ESSENZA basic software information	
Name and Version:	syngo MR E11
Date of release:	May 20, 2016
Size in GB:	5.56 GB

5.20. *MRI-Stretcher, non-magnetic:*

Basic MRI-Stretcher information	
Mass, kg	41kg
Dimensions, mm	1900x650x800mm (l x w x standard height)
Heights, mm	650,700,750 and 800
Diameter of wheels, mm	100
Height adjustment step	50 mm
Effort for trolley movement (no more than)	100 N
Effort for STOP ONN	No more than 150N
Removable elements	NA
Effort for wheels movement	0,35 N
Max load	No less than 160 kg
Gap between wheels and floor	No less than 3 mm

5.21. *PMU Wireless Physio Control; PMU Sensor Kit*

Basic PMU Wireless Physio Control; PMU Sensor Kit information	
Respiratory range	15-40 times/min
Pulse range	50-160 times/min

Frequency	2.4GHz
Modulation parameters	FM
IEC 60601-1 classification	Component with internal power supply Applied part: Type BF
Internal power supply	DC 3.7V
IP	IPX0
operating mode	Continues
MR compatibility	MR-safe
Respiratory cushion	
Manufacturer	Pi Products
PN	07389500
Weight	0.32kg
Dimensions	206mm × 64mm × 10mm
Application cushion	
Manufacturer	PolyForm GmbH
PN	10433187
Weight	0.4kg
Dimensions	425mm × 80mm × 20mm
Respiratory belt	
Manufacturer	Pi Products
PN	07389500
Weight	0.32kg
Dimensions	1390mm × 105mm
Pulse sensor	

Manufacturer	Mammendorfer Institut FÜR Physik UND Medizin GmbH
PN	10432918
Weight	0.6kg
Dimensions	58mm × 1380mm × 18mm
Charger station	
Manufacturer	Frör Kunststofftechnik GmbH
PN	10432957
Weight	0.604kg
Dimensions	140mm × 122mm × 33mm
Input supply	DC 5V
Input current	500mA
Power	2.5W
Batteries	
Weight	26 g
Dimensions (mm)	61 mm × 37.5 mm × 5.95 mm
Supply	DC 3.7V
Capacity	1320 mAh
Max. working time	No less than 12 hours
Max charging time	Approx 4 hours

5.22. Patient Data Display

Basic Patient Data Display information	
Manufacturer	Sasse
PN	10433484

Weight	0.9kg
Dimensions	350 mm × 220 mm × 35 mm

5.23. Keyboard

Basic Keyboard information	
Manufacturer	Fujitsu Technology Solutions
PN	10661060
Weight	2.25kg
Dimensions	465 mm×180 mm×20 mm
Cable length	4 m

5.24. Video Distributor for monitor

Basic Video Distributor for monitor information	
Manufacturer	SIEMENS
PN	10591047
Weight	5.3kg
Dimensions	340 mm×240 mm×50 mm

5.25. Trigger Converter

Basic Video Distributor for monitor information	
Manufacturer	SIEMENS
PN	10092360
Weight	0.7kg

Dimensions	7.681cm3
------------	----------

5.26. MR Workplace Table

Basic MR Workplace Table information	
Manufacturer	Dongguan weige steel supreme office equipment co., LTD
PN	10498411
Weight	92 kg
Dimensions	1020 mm x 910 mm x 720 mm

5.27. MR Workplace Container

Basic MR Workplace Container information	
Manufacturer	Dongguan weige steel supreme office equipment co., LTD
PN	10498412
Weight	40 kg
Dimensions	910 mm x 730 mm x 910 mm

5.28. MagnaCoil

Basic MR Workplace Container information	
Manufacturer	Magnacoustics Inc
PN	11234549
Weight	0.4kg
Dimensions	0.75" x 0.50" 0.75" x 0.38"

5.29. RF-Cabin

Basic RF-Cabin information

Manufacturer	SIEMENS
PN	03223448

5.30. RF-Traversal Beams

Basic RF-Traversal Beams information

Manufacturer	SIEMENS
PN	07014678 01141725
Weight	479.7kg 200kg
Dimensions	2.855m3 0.58 m3

5.31. Filter Cabinet, big

Basic Filter Cabinet, big information

Manufacturer	SIEMENS
PN	10131488 10498484 10838252 10837768
Weight	130 kg
Dimensions	1179 mm × 550 mm × 816 mm

5.32. Adapter, RF Cabin

Basic Adapter, RF Cabin information

Manufacturer	SIEMENS
--------------	---------

PN	10277265
Weight	0.35kg
Dimensions	3.808dm3

5.33. Water hose extension

Basic Water hose extension information	
Manufacturer	Dimplex (Shenyang) Electric Co., Ltd
PN	10352392
Weight	119.5kg
Dimensions	1012 mm ×680 mm ×940 mm

5.34. Seismic Kit

Basic Seismic Kit information	
Manufacturer	Shenzhen Deyihua Machinery Fittings Co. Ltd.
PN	10840554
Weight	41.5kg
Dimensions	0.192m3

5.35. UPS

Basic UPS information	
Manufacturer	Emerson
Model	GXT4
PN	11291557 UPS 1ph 3k VA GXT4 + power cable 11291559 Battery module UPS GXT4
Weight	40kg

Dimensions	408 mm × 430 mm × 85mm
------------	------------------------

5.36. *Distribution board*

Basic Distribution board information	
Manufacturer	Eaton
Model	USV-SYS-543-390-V1-0-A1
PN	08114717
Weight	5.0kg
Dimensions	510 mm × 330 mm × 320 mm

5.37. *ECG-electrodes*

Basic Distribution board information	
Manufacturer	ConMed
Model	CLEARTRACE ²
PN	07437861
Weight	0.02 kg
Dimensions	54mm × 55mm
Shelf life	2 years

5.38. *Service Syphon Kit*

Basic Service Syphon Kit information	
Manufacturer	Siemens HC Ltd. Magnet Technology
PN	10126651
Weight	5.15kg
Dimensions	0.076m3 (Volume)

5.39. HE-SYPHON DEWAR END

Basic HE-SYPHON DEWAR END information	
Manufacturer	Siemens HC Ltd. Magnet Technology
PN	01564207(HE-SYPHON DEWAR END 3.5M) 01194344 (HE-SYPHON DEWAR END 1.6M)
Weight	01564207: 15.3kg 01194344: 13kg
Dimensions	01564207: 0.306m3 (Volume) 01194344: 3.175Hectoliter(Volume)

5.40. SYPHON EXT.LEG

Basic SYPHON EXT.LEG information	
Manufacturer	Siemens HC Ltd. Magnet Technology
PN	01562870
Weight	0.08kg
Dimensions	1.134 dm3(Volume)

5.41. Laser positioning

Basic Laser positioning information	
Length wave	650 nm
relative intensity	$\leq 25 \text{ W/m}^2$
Max power	0.6 mW

5.42. LCD Color Remote Monitor:

Basic LCD Color Remote Monitor information	
--	--

Resolution, pixels	1280 x1024
Color numbers	16.7M
Brightness, cd/m2	Typical 300
Display size	19"
Manufacturer	Siemens A&D SE DT/ EIZO GmbH
PN	10683899
Weight	≤ 3.8 kg
Dimensions	277 mmx 294mm x 80mm
Supply	120/230 VAC +/- 10% automatic selection
Frequency	50/60 Hz

6. The list of components and accessories.

The magnetic-resonance system MAGNETOM ESSENZA with accessories		
##	Name of component	Description
BASIC BLOCK:		
1.	MAGNETOM ESSENZA system – the whole body magnet.	Whole-body superconductive 1.5T magnet with active shielding technology with counter coils.
2.	Tim [25x8] V-engine #Ez Tim [46x16] V-engine #Ez	The whole-body gradient system engine.
3.	MAGNETOM ESSENZA system – integrated Body Coil.	Integrated no tune transmit/receive Body Coil.
4.	MAGNETOM ESSENZA system – Front cover.	Magnet cover.
5.	Patient table w/o vert.move,8-ch #Ez Patient table w vert.move, 8-ch. #Ez Patient table w/o vert.move,16ch #Ez Patient table w vert.move,16-ch. #Ez	The patient table is mounted directly to the magnet assembly. The table comes with a comprehensive set of patient positioning cushions.
6.	MAGNETOM ESSENZA system – Electronic cabinets.	Electronics cabinets: - Gradient cabinet; - Control cabinet; - System cooling unit. Function: Interface between the on-site water chiller (of any brand/type) or Interface to the central hospital chilled water supply.
7.	Passive Cooling System,8-channel#Ez Active Cooling System, 8-channel#Ez Passive CoolingSystem,16-channel#Ez Active Cooling System,16-channel#Ez Fan Unit, Standard #Ez	Water cooling solution.

8.	MAGNETOM ESSENZA system – the alarm box.	The alarm box has the following functions: - Displays alarm signals; - Switches the MR system on and off; - Magnet Stop/ Magnet Quench.
9.	MAGNETOM ESSENZA system – emergency switches.	The MR system has different types of emergency switches. (E.g. in case of accidents with attracted metal parts, fire, in case of accidents or injury due to table movement).
10.	MAGNETOM ESSENZA system – the intercom.	The intercom allows personnel and patients to communicate during the examination.
11.	syngo Acquisition Workplace.	syngo Acquisition Workplace is the workplace in the control room. It includes the host processor with the operating elements monitor and mouse.
12.	PC Keyboard US English #Tim PC Keyboard US english #Ez	The MR system is equipped with an original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.
13.	The application package “Tim Application Suite”: -Basic software license: basic license key for running the syngo software; -Basic MR applications license; -3D functionality; -Argus (cardiac) basic viewing; -DICOM worklist support; -PhoenixZIP; -Different standard measurement/examination and post-processing functionality.	Provides a complete range of clinically optimized examinations for all regions. Enables easy reproduction of measured programs and therefore repeated examinations are efficiently facilitated.
14.	Head Matrix Coil #Ez	The Head Matrix Coil is a fully iPAT-compatible no-tune coil. The Head Matrix Coil can be used e.g. for head examinations, MR Angiography, combined head/neck examinations (in combination with the Neck Matrix Coil) or for imaging of the TMJ (temporo mandibular joints).
15.	Neck Matrix Coil #Ez	The Neck Matrix Coil is a fully iPAT-compatible no-tune coil. The upper coil part can be removed for better patient handling. The lower coil part can remain on the

		table for most examinations. The Neck Matrix plugs into the Head Matrix for combined head-neck imaging. Together with the Isocenter Spine Matrix it can be used for applications like neck or cervical spine examinations, imaging of the Larynx/Esophagus and Mediastinum MR Angiography, and combined head/neck examinations. It thus works as a neurovascular coil.
16.	MAGNETOM ESSENZA system – the isocenter spine matrix coil.	Isocenter Spine Matrix is a fully iPAT-compatible no-tune coil which is located under the patient table in the isocenter of the magnet. The unique concept is designed to optimize system utilization and workflow.
17.	The Flex coil interface.	The Flex Coil Interface is required for connecting the Flex coils, the Loop coils and the Endorectal coil.
18.	MAGNETOM ESSENZA system – flexible multi-purpose coils (2 pcs.).	Light-weight, very flexible, iPAT compatible, no-tune receiver coil made of soft, smooth material. It is used for examinations of the upper and lower extremities (e.g. small to medium shoulder, wrist, elbow or ankle).
19.	MAGNETOM ESSENZA system – connection cables.	Connection cables.
20.	MAGNETOM ESSENZA system – quality calibration phantoms.	Measurement phantoms are used for quality measurements.
21.	Phantom holders.	Phantom holders for transportation and storage.
22.	Russian User Documentation.	Instruction for use.
23.	MRI-Stretcher, non-magnetic	MRI-Stretcher, slightly magnetic, chrome-nickel compounded for patient transport inside the clinic building. Continuously adjustable head piece with spring brake, synthetic leather cushion, movable barrier on either side, stretcher movable on four fixable wheels diameter 100 mm.
24.	IV Stand, non-magnetic	IV-Rod, fix mounted, for AGA non magnetic stretcher. - Diameter: 180 mm; - Length: 1720 mm; - With four hooks.
25.	Glider Positioning Kit	Positioning System for patient transport between bed or wheelchair and MR-System. These sliding boards consist of a special material with very low coefficients of friction. They can be folded, take up little space and can be transported easily. They are firm but flexible and adapt to different shapes of the body.
26.	Instrument table, non-magnetic	Mobile 4-leg frame construction made of steel. Two shelves: surface ground to reduce glare, resistant to disinfectants, removable. 4 plastic swivel castors, two of

		them with brakes, wear resistant, non-discoloring when in contact with other material, electrically conductive.
27.	ECG-electrods	ECG-electrods
28.	Flex Large	Light-weight, very flexible, iPAT compatible, no-tune receiver coil made of soft and smooth material. It is used for examinations of extremities (e.g. medium to large shoulder, hip or knee) and for abdominal examinations.
29.	Flex Small	Light-weight, very flexible, iPAT compatible, no-tune receiver coil made of soft and smooth material. It is used for examinations of the upper and lower extremities (e.g. small to medium shoulder, wrist, elbow or ankle).
30.	PMU Wireless Physio Control #Ez PMU Sensor Kit #Ez	Physiological Measurement Unit (PMU) - Wireless Physio Control for wireless triggering, synchronizes the measurement with the physiological cycles of cardiac and/or respiratory motion.

ACCESSORIES:

1.	syngo MR Workplace.	The syngo MR Workplace allows for evaluation, documentation and post-processing of previously measured images while acquiring images at the syngo Acquisition Workplace.
2.	The additional workstation based on User solution.	The additional workstation based on User solution entry level sales package for advanced visualization customers. It provides the same software functionality, as the conventional syngo.via product.
3.	In-Room syngo Acquisition Workplace	The In-Room syngo Acquisition Workplace (In-Room syngo Acq WP) is an additional operating console in the examination room. It is used for image viewing and MR system operation.
4.	Optional Dot Engine Workflows.	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.
5.	Optional Application Packages	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.
6.	Optional Postprocessing Packages.	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.

7.	Additional optional single Applications	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.
8.	Argus Flow; Argus Dynamic Signal	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.
9.	syngo Expert-i (syngo Expert-i MRWP #T+D)	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.
10.	Syngo Security package	Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. This package supports customers in their achieving compliance with the HIPAA "Privacy" rule.
11.	16-channel Head/Neck Coil #Ez	The 16-channel Head/Neck Coil is a fully iPAT-compatible no-tune coil. Together with the Isocenter Matrix it can be used for applications like neck or cervical spine examinations, imaging of the Mediastinum, MR Angiography and combined head/neck examinations. It thus works as a neurovascular coil.
12.	Body Matrix Coil #Ez	The multi-element Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio. Applications: - Thorax (incl. heart); - Abdomen; - Pelvis; - Hip.
13.	PA Matrix Coil #Ez	The multi-element Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio. Applications: - High-resolution angiography of both legs incl. pelvis with highest signal-to-noise ratio; - Visualization of the iliac arteries and aorta.
14.	Extremity Matrix Coil #Ez	Extremity Matrix Coil is an iPAT compatible no tune receive coil for knee joint examinations. The anatomically adapted coil allows scanning with comfortable

		positioning of both legs. The upper part of the coil can be removed for easy patient positioning.
15.	14-channel Extremity Coil #Ez	Novel 14-channel high-performance extremity coil, optimized for knee imaging. Main features: - 18-element design (3 rings with 6 coil elements in head-feed direction) with 14 integrated preamplifiers; - iPAT-capable in all directions. Thanks to its 14-channel design this coil is perfectly suited for high-resolution images with excellent SNR. With the arrangement of the antennas in three rings of 6 elements each, the coil is specially designed for parallel imaging with high acceleration factors in all directions. The coil is positioned on a laterally movable support and therefore allows for comfortable patient positioning of both legs for off-center examinations. Furthermore, the upper part can be removed for easier patient positioning. Additional cushions allow for optimum patient immobilization. The integrated hybrid transmission function makes single knee excitation with reduced RF power possible. Besides that, it eliminates wraparound artifacts from the contralateral knee and flow artifacts from upstream vessels.
16.	8-channel Foot/ Ankle Coil #Ez	The 8-channel foot-ankle coil is an iPAT-compatible "no-tune" receiver coil for the examination of the foot and the ankle joint.
17.	8-channel Wrist Coil #Ez	The 8-channel wrist coil is an iPAT-compatible "no-tune" receiver coil for the examination of the wrist and fingers.
18.	Focus Shoulder Array Coil #Ez Focus Shoulder Arr.Coil, small #Ez Focus Shoulder Coil package #Ez	The Focus Shoulder Array Coil is an iPAT compatible coil for examinations of the left or right shoulder. The Focus Shoulder Array Coil features: - No coil tuning - Focus Wing with integrated shim device focuses the homogeneity volume of the magnet to the off-center anatomy. - For narrow or wide shoulders the coil can be attached at different positions on the base plate. - One base plate pad and one head rest for high patient comfort are included.
19.	Breast Matrix Coil #Ez	The breast matrix coil delivers brilliant image quality for high-resolution 2D and 3D MR breast imaging.
20.	Endorectal Coil Interface #Ez	The endorectal coil is used for prostate, cervix or colon MRI.

21.	4-channel Flex Coil Interface #Ez	Interface with 4 integrated preamplifiers for the connection of the following flexible coils to the system: - 4-channel Flex Coil, large; - 4 channel Flex Coil, small; - Endorectal Coil.
22.	Special Purpose Coil #Ez	The special purpose coil is a no-tune receiver coil for the examination with small Field-of-Views.
23.	Cushions and positioning aids	Cushions and positioning aids
24.	Set of positioning aids.	Set of positioning aids
25.	Head Matrix accessories Coil #Ez	Head Matrix Coil accessories for patient positioning
26.	Neck Matrix Coil accessories #Ez	Neck Matrix Coil accessories for patient positioning
27.	16-channel Head/Neck Coil accessories #Ez	Head/Neck Coil accessories for patient positioning
28.	Body Matrix Coil accessories #Ez	Body Matrix Coil accessories for patient positioning
29.	Breast Matrix Coil accessories #Ez	Breast Matrix Coil accessories for patient positioning
30.	PA Matrix Coil accessories #Ez	PA Matrix Coil accessories for patient positioning
31.	Focus Shoulder Array Coil accessories #Ez Focus Shoulder Arr.Coil accessories, small #Ez Focus Shoulder Coil accessories package #Ez	Shoulder Coil accessories for patient positioning

32.	Extremity Matrix Coil accessories #Ez	Extremity Matrix Coil accessories for patient positioning
33.	14-channel Extremity Coil accessories #Ez	14-channel Extremity Coil accessories for patient positioning
34.	4-channel Flex Coil, large accessories #Ez 4-channel Flex Coil, small accessories #Ez	Flex Coil accessories for patient positioning
35.	Endorectal Coil Interface accessories #Ez	Endorectal Coil accessories for patient positioning
36.	8-channel Foot/ Ankle Coil accessories #Ez	Foot/ Ankle Coil accessories for patient positioning
37.	8-channel Wrist Coil accessories #Ez	Wrist Coil accessories for patient positioning
38.	Coil Storage Cart #T+D,Ez PA Matrix Coil Storage Cart #Ez	Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories. The cart may be rolled to convenient locations in the examination room and can be opened up to work like a shelf.
39.	Patient Supervision TV #Ez	This package contains a special video camera for monitoring the patient during an MR examination, conveniently mounted on the wall of the examination room. The information is displayed at an LCD monitor in the control room.
40.	LCD Color Remote Monitor	This remote monitor is used for connection in parallel to the syngo Acquisition Workplace monitor. The LCD remote monitor is of the same type as the LCD color monitor of the syngo Acquisition Workplace with the same high-quality characteristics: - can be tilted horizontally forward and backward; - integrated gamma correction for optimum display of radiographic grayscale; - automatic backlight control for long-term brightness stability.

41.	Patient Data Display #Ez	LCD monitor mounted on the magnet cover for display of physiologic (ECG, pulse, respiration) signals in the examination room. Including front cover with outlet for the LCD display.
42.	Color monitor	Color monitor
43.	PC Keyboard US English #Tim	Additional PC keyboard.
44.	Video Distributor for monitor #Ez	Video distributor for connection of up to two additional color LCD monitors as well as an In-Room Acq WP syngo (in the examination room parallel to the syngo Acquisition Workplace operator monitor).
45.	Trigger Converter	An optical trigger signal is available to trigger external stimulation devices in f experiments.
46.	MR Workplace Table #Ez,Am,Se	Ergonomically designed table for: - Color monitor; - Keyboard; - Mouse; - Patient communication unit; - Patient supervision display.
47.	MR Workplace Container #Ez,Am,Se	Extra case for the syngo host computer with sliding front door to allow change of storage media.
48.	MagnaCoil	The MagnaCoil Headset can be used as an alternative headphone (for example, in head coils) to the standard Siemens headphones.
49.	RF-Cabin #1.5T	Steel-full-room RF cabin, in a module method of building. Standard outside dimensions: L x W x H =6,590 x 4,090 x 2,750 mm
50.	RF-Traverse Beams	4 traverse beams for the RF cabin.
51.	RF-Interior, LED	Standard interior decor for the RF cabin.

52.	RF Cabin Sound Proofing	Sound protection system. Due to building requirements according to DIN 4109 (sound protection for structural engineering) it is possible to get the following values: - Sound protection, walls > 52 dB(A); - Sound protection, ceiling > 56 dB(A); - Sound protection, observation window > 50 dB(A).
53.	Filter Cabinet, big	Cabinet for filter plate inside of the RF cabin. Wall-mounted cabinet custom-made. Outside dimensions appr. 1.300 x 420 x 2.700 mm (L x W x H).
54.	Threshold-free door	Soundproof threshold-free door for the RF-Cabins that can be opened and closed without effort by pneumatic mechanism.
55.	Emergency hatch	With decorative plate on both sides, white, to be opened from the outside.
56.	Observation/ Day-light window	Observation / daylight window that can be opened, installation sideways.
57.	Adapter, RF Cabin	2 adapter plates for the installation of the (smaller) filter plate into the (larger) filter plate opening of an existing RF cabin.
58.	Water hose extension 20m #Ez	20 m water hose extension including connectors.
59.	Seismic Kit	It comprises: - 3/4" tufnol spacers which have hole offsets from zero to 4.2 mm in 0.7 mm increments and washers; - Radial gap checking tool; - Vertical gap setting tool; - 1" tufnol spacer kit - only to be used with 1" floor anchors; - Loctite 243 (10 ml bottle) - 4 x M20 washers.
60.	UPS	Uninterruptible Power Supply

61.	Distribution board	Distribution board
62.	Metal detector	Metal detector
63.	A fire extinguisher	A fire extinguisher
64.	Service Syphon Kit	Service Syphon Kit
65.	HE-SYPHON DEWAR END	HE-SYPHON DEWAR END
66.	SYPHON EXT.LEG	SYPHON EXT.LEG
67.	Floor plate for table	Base Plate for the fixed table
68.	Adhes. Plastic-Mastic complete	Adhesive Plastic-Mastic complete kit consisting of adhesive component A, adhesive component B, gloves and spatula.
69.	Respiratory cushion	Respiratory cushion with pressure hose
70.	Application cushion	Application cushion
71.	Respiratory belt	Respiratory belt
72.	Charging station	The charging station is installed separately near the syngo Acquisition Workplace and is used for storing both units.

7. Capabilities of the Medical Device and ways of integration with other medical devices.

Name of the device	Manufacturer
MRI-Stretcher, non-magnetic	AGA
IV-Rod, non-magnetic f. Strecher	AGA
Instrument table, antim agnetic	AGA
MRT Pat Stretcher, hight adjustable	AGA

MR Injector Spectris Solaris EP (D)	BAYER MEDICAL CARE
MR Injector Spectris Solaris EP	BAYER MEDICAL CARE
MRXperion Injector (DE)	BAYER MEDICAL CARE
MRXperion Injector (RoW)	BAYER MEDICAL CARE
EMPOWER MR Injector	BRACCO
MR-ECG Electr.Conm ed 2700/600pcs .	CONMED
Draeger Fabius MRI, CE	DRAEGER
Draeger Fabius MRI, CE, Pin Index	DRAEGER
DORADOnova3 MR3 Bridge (red)	LAP
DORADOnova3 MR3 Bridge (green)	LAP
Laser Bridge LAP DORADOnova3 red	LAP
Laser Bridge LAP DORADOnova3 green	LAP
MagnaPlug Standard (250 pair)	Magnacoustics
MagnaPlug Standard (500 pair)	Magnacoustics
MagnaPlug Mini (250 pair)	Magnacoustics
MagnaPlug Mini (500 pair)	Magnacoustics
Medtron Accutron MR	MEDTRON
Medtron Accutron MR (EU)	MEDTRON
Accutron MR3 Injector (DE)	MEDTRON
Accutron MR3 Injector (RoW)	MEDTRON
MRI-Wheelchair, non-magnetic	MEYRA
MIPM TeslaM3 Basic	MIPM
MIPM TeslaM3 Capno	MIPM
MIPM TeslaM3 Basic Anesthesia	MIPM
TeslaDUO Configuration: SpO2 & NIBP	MIPM
TeslaDUO Configuration: SpO2 only	MIPM
TeslaDUO Configuration: NIBP only	MIPM
TeslaDUO IR Remote Control for NIBP	MIPM
TeslaDUO Trolley	MIPM
NNL Basic fMRI System+40" Monitor	NNL NORDICNEUROLAB
NNL MR In-Room-Monitor 40"	NNL NORDICNEUROLAB
NNL MRI Imaging In-Room-Monitor 40"	NNL NORDICNEUROLAB
Variety PositioningAids Orthopedics	Noras
Variety Positioning Aid Pediatrics	Noras
Ohropax Soft (100 pairs)	OHROPAX
Ohropax Classic (100 pairs)	OHROPAX
IV Stand, non-magnetic	PROVITA
HW starter set for Elastography	RESOUNDANT
MAGLIFE light S/N	SCHILLER
MAGLIFE Serenity E/S/N - SpO2	SCHILLER
MAGLIFE Serenity E/S/N - SpO2 WLAN	SCHILLER
IBP Invasive Blood Pressure, 2 ch.	SCHILLER
Temperature Option	SCHILLER
Spacelabs MR anesthesia w. suction	SPACELABS

The Glider Positioning Kit	SWERECO
Fire extinguisher, anti-magn., emp.	TOTAL
Fire Extinguisher, non-magnetic	TOTAL
Foot Step, non-magnetic	TRUMPF
Ulrich MR-Injector, mississippi	ULRICH
Ulrich MR-Injector tennessee	ULRICH
Ulrich MR-Injector, mississippi (D)	ULRICH
Ulrich MR-Injector tennessee (D)	ULRICH
SpaceStation / Perfusor Space	B.Braun
Earplugs	Bilsom / Honeywell
Sphygmomanometer	Erka
Non- Magnetic Stethoscope	Erka
MRI compatible cleaning containers	Hara
MRI Anesthesia Cart	Harloff
MRI Anaesthesia System Leoni	Heinen & Loewenstein
Esys Patient Entertainment+ fMRI System	Invivo
MRI Infusion Pump	Iradimed
Magnetic Laryngoscope	MRI MED
Earplugs „MiniMuffs“	Natus Medical
Communication System	Optoacoustics
Compatible Containers	Rainer Medizintechnik
MRI Prism Glasses	ScanSound
Chair	Score

8. Labeling information.

8.1. At electronic cabinets.

Name plate label:

MAGNETOM ESSENZA

(01) 04056869006772 GTIN
(240) 10837633 MODEL
(422) 276 ORIGIN DE
(21) XXXXXX SERIAL

SYSTEM IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in Germany

Model	Model
SERIAL	Serial number
Made in Germany	Made in Germany
SYSTEM	System
	Not to throw out. To give to utilization point.
	Manufacturing date

Revision label:

Revision:															Model - No.:																	
00	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65
66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98

Patent label:

This product is covered by one or more of the following US-patents:			
4892788	8287837	7038463	7477067
5034882	8300781	7060817	7482158
5138268	8307374	7061288	7485437
5189370	8335820	7063744	7485444
5230090	8338332	7084828	7501828
5284888	8338717	7078901	7511488
5304828	8361123	7081720	7528804
5307014	8388080	7082834	7541904
5308107	8388588	7088868	7580830
5328288	8386480	7088718	7570058
5337000	8400163	7108716	7673288
5461877	8404185	7118543	7674371
5468401	8414488	7128334	7678837
5487017	8420870	7146338	7688737
5482124	8433548	7164270	7809058
5512828	8434412	7165297	7809058
5523888	8482547	7184288	7822823
5541614	8501273	7184333	7833281
5548081	8515478	7181848	7833282
5560471	8625537	7183865	7834827
5574372	8638248	7187172	7839010
5581184	8646478	7183418	7848363
5614827	8658855	7206783	7872488
5623207	8658841	7218383	7889282
5633588	8673718	7218108	7889270
5682112	8689431	7218114	7884386
5688474	8681148	7224185	7702376
5708381	8834831	7228258	7702378
5777475	8838038	7230427	7704650
5860143	8840018	7233148	7728588
5886624	8860925	7235874	7728891
5892838	8887817	7238141	7737787
5898688	8889068	7239144	7747313
5900825	8708373	7242401	7754875
6023788	8711434	7246125	7791340
6026718	8711738	7246762	7808428
6043861	8714083	7248047	7818818
6084204	8717487	7260785	7821288
6084408	8734872	7263818	7840048
6084410	8782805	7263828	7847554
6101844	8771872	7264435	7852829
6107788	8781378	7278882	7859281
6111412	8784872	7286008	7874031
6114854	8786037	7288823	7884808
6133737	8801835	7318234	7884807
6146028	8822445	7318324	7885442
6164088	8825885	7331888	7888884
6167181	8838722	7333848	7891188
6167183	8841888	7340225	7894782
6180400	8841888	7342388	7824006
6180463	8850208	7386538	7824807
6188408	8864883	7387838	7827442
6170018	8880977	7388915	7848541
6186031	8882161	7372288	7885082
6188287	8818721	7388875	7888066
6206348	8830482	7387248	7870888
6216811	8837018	7417427	7877846
6218338	8843551	7421103	7886142
6238204	8852087	7428882	8000770
6238208	8877430	7432713	8008818
6240310	8882588	7432784	8018403
6262401	8882485	7438742	8030827
6263101	7002344	7440781	8055188
6265872	7008388	7468848	8062728
6281765	7038814	7488860	8083878
6288647	7034531	7474877	8086045
6287833	7038451	7474812	8138758

Artwork of the MD labeling in Russian language*

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~3/PE 380/400/420/440/460/480 В

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*This label isn't applicable for the systems which are delivered before Feb 2020.

Gradient system V-engine:

Grad_coil_011_kpl

(01) 04056869225081	GTIN
(240) 10164200	MODEL
(422) 156	ORIGIN CN
(21) 4344	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaixin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

SIEMENS

GPA 16_13_D

(01) 04056869092133	GTIN
(240) 11106650	MODEL
(422) 276	ORIGIN DE
(21) 1717	SERIAL

IVK


Siemens Healthcare GmbH, Henkestr. 127
91052 Erlangen, Germany

Made in Germany

GTIN	Global product number GTIN
Model	Model
ORIGIN DE	Country of Origin (Germany)
SERIAL	Serial number
Made in Germany	Made in Germany
Made in China	Made in China
	Manufacturer

Patient table


Patient table without vertical move,8-ch

Patient table with vertical move, 8-ch.

Patient table without vertical move, 16ch

Patient table with vertical move,16-ch

GTIN	Global Product number GTIN
Model	Model
ORIGIN DE	Country of Origin (Germany)
SERIAL	Serial number
Made in China	Made in P.R.China
	Manufacturer

Magnet cooling system

 GLEN DIMPLEX DEUTSCHLAND		designed by Glen Dimplex Deutschland GmbH Am Goldenen Feld 18 D-95326 Kulmbach Germany	SDE Made in China	
Model-No. / no. modèle	2NK6 770 K			
consisting of. / consisten en	2NK6 770 W + 2NK6 770 A			
Voltage / tension	460 V			
Frequency / fréquence	60 Hz			
Phase / phase	3			
Wire / fil	3			
Max Fuse Size / capacité fusible	50 A			
Min Circuit Ampacity / intensité min. admissible	36 A			
Ampere Primary Pump / ampère primaire pompe	4.0 A			
Ampere Secondary Pump / ampère secondaire pompe	4.0 A			
Ampere Outdoor Fan / ampère ventilateur	2 x 1.1 A			
Ampere Compressor / ampère compresseur	12 A			
Ampere Cryogenic Unit / ampère cryogénique unité	13 A			
Refrigerant / frigorigène	134a			
Amount / quantité	6.4 lb			
High Side Pressure / zone haute pression	294 psig			
Low Side Pressure / zone basse pression	220 psig			
Date Code / code fabrication	FD			
 PROCESSING WATER COOLER				
REV: 00 01 02 03 04 05 06 07 08 09 10				454507.22.02

Name plate Active Cooling System

 GLEN DIMPLEX DEUTSCHLAND		designed by Glen Dimplex Deutschland GmbH Am Goldenen Feld 18 D-95326 Kulmbach Germany	SDE Made in China	
Model-No. / no. modèle	2NK6 771			
Voltage / tension	480 V			
Frequency / fréquence	60 Hz			
Phase / phase	3			
Wire / fil	3			
Max Fuse Size / capacité fusible	30 A			
Min Circuit Ampacity / intensité min. admissible	19 A			
Ampere Pump / ampère pompe	4.0 A			
Ampere Cryogenic Unit / ampère cryogénique unité	13 A			
Medium / liquide Type / type	water / eau			
Cooling water High Side Pressure / eau de refroidissement zone haute pression	88 psig			
Fresh water High Side Pressure / eau douce zone haute pression	88 psig			
Date Code / code fabrication	9905			
 CULUS		WATER CHILLER		
REV: 00 01 02 03 04 05 06 07 08 09 10				454508.22.02

Name plate Passive Cooling System

GTIN	Global product number GTIN
Model	Model
ORIGIN DE	Country of Origin
SERIAL	Serial number
Made in China	Made in China

	Manufacturer
	Not to throw out. To give to utilization point.
	EC Conformity

MAGNETOM ESSENZA system – the intercom


Model	Model
SERIAL	Serial number
Made in Germany	Made in Germany
SYSTEM	System
	EC Conformity


Head Matrix Coil

Neck Matrix Coil

MAGNETOM ESSENZA system – the isocenter spine matrix coil


The Flex coil interface


16-channel Head/Neck Coil

Body Matrix Coil

PA Matrix Coil

Extremity Matrix Coil


14-channel Extremity Coil

8-channel Foot/ Ankle Coil


8-channel Wrist Coil


Focus Shoulder Array Coil

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
Breast Matrix Coil


GTIN	Global product number GTIN
Model	Model
ORIGIN CN	Country of Origin (China)
ORIGIN DE	Country of Origin (Germany)
SERIAL	Serial number
Made in Germany	Made in Germany
Made in China	Made in China
SYSTEM	System

	Manufacturer
---	--------------

Endorectal Coil

Endo Interface

Endo adapter

4-channel Flex Coil

4-channel Flex Coil, small

4-channel Flex Coil, large


Special Purpose Coil


GTIN	Global product number GTIN
Model	Model
ORIGIN CN	Country of Origin (China)
SERIAL	Serial
Made in Germany	Made in Germany
Made in China	Made in China
Made in USA	Made in USA

	Manufacturer
	EC Confirmation

PMU Wireless Physio Control

PERU ECG Respiratory Unit

PPU Pulse Unit

Model	Model
SERIAL	Serial number
	Applied part type BF

Trigger Converter

GTIN	Global product number GTIN
Model	Model
ORIGIN DE	Country of Origin (Germany)
SERIAL	Serial number
Made in Germany	Made in Germany
	Manufacturer

UPS

Model	Model
SERIAL	Serial number
	Not to throw out. To give to utilization point.
	CE conformity

ECG electrodes

	Manufacturer
Made in USA	Made in USA
	Catalogue number

	Temperature range
	Single – use only
	LOT
	Attention! Please follow instruction for use
	Please follow instruction for use
	EC confirmation
	User before

HW

Monitor

GTIN	Global product number GTIN
Model	Model
ORIGIN CN	Country of Origin (China)
ORIGIN CZ	Country of Origin (Czech Republic)
SERIAL	Serial number
Made in China	Made in China

Keyboard

GTIN	Global product number GTIN
Model	Model
ORIGIN CZ	Country of Origin (Czech Republic)
SERIAL	Serial number
Made in China	Made in China
SYSTEM	System
	Not to throw out. To give to utilization point.
	EC Confirmation

For laser light localizer:

Safety certificate label (international)

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The EAC label:

The Eurasian Conformity mark is a product marking to indicate that the product conforms to all appropriate Technical Regulations of the Customs Union assessment procedures.

8.2. Label information:

SIGNS AND SYMBOLS		
		Magnet Stop Press this button to remove the magnet field. Emergency use only!
		Table Stop This symbol indicates the Table Stop button. It is located on the intercom, on the patient table and optionally on the magnet cover.
	Attention, magnetic part, observe operator manual This label is located on the cover of the MR system (bottom left).	
	Trigger connection This label is located on the cover of the MR system (bottom left).	

	Maximum patient weight This label is located on the patient table.
---	--

SIGNS AND SYMBOLS	
	Maximum total weight (patient + patient table) This label is located on the dockable patient table.
	Avoid cable loops, observe operator manual This label is located on the Endorectal Interface (3T).
	Start and Stop (the MR measurement) These labels can be affixed to the footswitches of the In-Room <i>syngo</i> Application Workplace.
	General warning sign Caution Attention: Consult accompanying documents. Indicates potential injury or damage to equipment. Observe safety precautions. ISO 7010-W001
	Caution: Strong magnetic field ISO 7010-W005

Consult instructions for use.

Follow operating instructions.

ISO 7010-M002

8.3. **Additional system characteristics:**

Magnet characteristics	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	< 0.1 ppm/h
Magnet bore diameter	82 cm*
Refill interval (typical)	10 years**
Boil-off rate (typical)	0.0 l/year**
Helium capacity	approx. 1400 liters
Cryostat	Stainless steel
Minimum helium level	40%

* without shim coils, gradient coil, RF body coil

**for typical clinical use, depending on sequences and operating time

RF system	
RF transmit coils	Body coil
Peak power of transmitter amplifier	15 kW
Transmitter bandwidth	500 kHz
Receiver bandwidth	500 Hz– 1 MHz (for each channel)
Max. applied RF field B1+	30 μ T
Max. specified B1+rms	6.7 μ T
Distance to isocenter, where RF transmit field is reduced by	
3 dB	0.18 m
10 dB	0.62 m

Gradient system (Tim 25 × 8)(Tim 46 × 16) V-engine	
Performance	
Max. amplitude	30 mT/m
Min. rise time	300 µs
Max. slew rate	100 T/m/s
Vector gradient performance (vector summation of all 3 gradient axes)	
Max. eff. amplitude	52 mT/m
Max. eff. slew rate	173 T/m/s

Patient table	
Width	55.5 cm (top)
Depth	210.5 cm
Max. patient weight	200 kg
Vertical table movement*	
Range	55 - 89 cm
Max. vertical movement time	13 s
Longitudinal table movement	
Max. range	205.5 cm
Max. speed	20 cm/s
Positioning accuracy	± 0.8 mm

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9. The list of International Standards and Directives applied to the Medical Device.

IEC60601-1:2005+CORR.1(2006)+CORR.2(2007+AM1(2012)	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2:2007	Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests
IEC 60601-1-6:2010	Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability
IEC 60601-2-33:2010/A2:2015	Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
IEC 62304:2006/ AMD 1 2015	Medical device software - Software life-cycle processes
IEC 62366:2007+A1:2014	Medical devices - Application of usability engineering to medical devices
ISO 14971:2007	Medical devices — Application of risk management to medical devices
ISO 13485:2016	Medical devices Quality management systems Requirements for regulatory purposes
ISO 10993-1:2009	Biological evaluation of medical devices: Part 1 Guidance on selection of tests
ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

10. Sterility.

The Medical device is non-sterile.

11. Disposal.

On disposing of the system or parts thereof, currently valid environmental legislation must be observed. Examples of environmentally relevant components are:

- Accumulators and batteries
- Transformers
- Capacitors
- Monitor picture tubes
- Phantoms

For details contact your local customer service representative or your Siemens regional office.

NOTICE: *System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.*

After usage the product can be classified as a medical waste of class "A".

All medical waste shall be collected, used, decontaminated, placed, stored, transported, listed and disposed of in accordance with the laws on public health and epidemiological welfare.

12. Transportation requirements and conditions.

Transport specifications - magnet	
Ambient temperature	-25°C to +60°C
Relative humidity	10% to 95% non-condensing
Atmospheric pressure	20 kPa to 106 kPa
Transport specifications for the system components	
Ambient temperature	-25°C to +60°C
Relative humidity	10% to 95%
Atmospheric pressure	70 kPa to 105 kPa

13. Storage requirements and conditions.

Storage specifications - magnet	
Ambient temperature	-25°C to +60°C
Relative humidity	10% to 95% non-condensing
Atmospheric pressure	20 kPa to 106 kPa
Storage specifications for the system components	
Ambient temperature	-25°C to +60°C
Relative humidity	10% to 95%
Atmospheric pressure	70 Pa to 105 kPa

14. Additional Maintenance information.

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Call center (24 hours) Customer's support center:

Tel: 8 800 550 50 22

Email: service.ru@siemens-healthineers.com

15. Cleaning and disinfection;

All instructions in the operator manual regarding cleaning, and when applicable regarding disinfecting and sterilization must be always observed.

Cleaning the LCD monitor/video display

The LCD monitor of the syngo Acquisition Workplace and the video display are cleaned in the same way.

- 1 Clean the LCD monitor/video display every other month.
- 2 Clean the LCD monitor/video display using a microfiber cloth.
- 3 In case the LCD monitor/video display are quite contaminated: Use a window cleaner.
- 4 Immediately remove water drops on the LCD monitor/video display.
- 5 Avoid scratching the surface area of the LCD monitor/video display.
- 6 Avoid impacts to the LCD monitor/video display.

The LCD monitor/video display are highly sensitive to mechanical damage.

Cleaning the video camera

- 1 Clean the video camera housing with a clean damp cloth.

- 2 Clean the optics of the video camera using a special lens cleaner.

Cleaning the data carrier

- 1 Clean contaminated data carriers with a clean cloth (cotton or microfiber).
- 2 Follow the manufacturer's notes when cleaning CD/DVD data carriers.

Cleaning the plugs and connectors

- 1 Dampen a soft cloth with water or a diluted household cleaner solution. Do not use organic solvents, such as alcohol or acetone.
- 2 Carefully wipe the plugs and connectors with the cloth. Do not touch the contacts.

Cleaning RF coils and positioning aids

- 1 Clean the RF coils and positioning aids using commercially available cleaning and disinfection solutions. Follow manufacturer's instructions.
- 2 Do not use hard or sharp objects (e.g. knives or tweezers) to remove residue.

Cleaning the patient table

- ◆ Clean the patient table using a liquid household agent. Use commercial disinfecting materials. However, do not use solutions with alcohol or acetone.

Cleaning the belt/straps of the patient table

- ◆ Wash the straps of the patient table at a temperature of 60 °C.

Cleaning and disinfecting the receptors

- 1 Do not use cleaners or disinfectants containing alcohol or ether.
- 2 Do not use hard or sharp objects (e.g. knives or tweezers) to remove residue.
- 3 Clean the receptors with a commercially available cleaner. Follow manufacturer's instructions.
- 4 Use a dampened cloth for cleaning.

Do not submerge the receptors in cleaning liquid.

- 5 Disinfect the receptors with a commercially available disinfecting agent. Follow manufacturer's instructions.

Disinfecting system components

Disinfecting sprays damage electronic components. For this reason, components may be cleaned with a semi-dry cloth only.

WARNING

Flammable cleaning or disinfection agents may cause fire or explosion!

Injury to patient

- ◆ Disinfect the system components with commercial disinfecting materials. However, do not use solutions with alcohol or acetone. Follow manufacturer's instructions.

Siemens does not test alcohol-containing, phenol-alkaline and alkaline-based disinfectants for harmful effects on surfaces. Avoid disinfecting system components with these materials.

Care and cleaning of floors

Do not use the following cleaning or care products:

- Sprays
- Silicon-based cleaning or care products
- Cleaning or care products with substances that release ammonia
- Cleaning or care products that destroy the anti-static properties of the floor covering
- ◆ Use commercially available cleaning or care products for the floor. Follow manufacturer's instructions.

16.Service life.

An expected system service life is **20 (twenty) years** in case of observance of the operation instruction and timely maintenance.

17.Warranty.

The warranty period supposed to be no less than 12 months from the date of commissioning but defined by the contract terms for every purchase.

18. Contacts for reclamation.

In case of any questions or problems with handling the system, please use the following contacts:

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC
Dubinskaya street 96
115093 Moscow,
Russia
Tel.: +7 495 737 10 00
Fax: +7 495 737 13 20
www.siemens-healthineers.com/ru/

Contacts for questions and reclamations in the country of Legal manufacturer:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center,
Gaoxin C. Ave. 2nd, Hi-Tech Industrial Park,
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86 0755 26525421
www.siemens.com.cn/ssmr

Glossary:

Name of the components according to Application for registration	Name of the components according to the technical and user documentation
Operator station: - control processor; - basic software; - monitor; - computer mouse.	syngo Acquisition Workplace
Head coil (Head Matrix)	Head Matrix
Neck coil (Neck Matrix)	Neck Matrix
Spine coil (IsoCenter Matrix).	IsoCenter matrix coil or ISO Matrix.
Flex multifunctional coils (Flex Small/Large)	Flex Small/Large
System of documentation and processing of received images, additional: - HW; - computer mouse; - color monitor.	syngo MR Workplace

Operator station of scanning and processing of the images in scanning room (HW, software, monitor, keyboard, mouse).	In-Room syngo Acquisition Workplace
Coil for Head/Neck (Head/Neck 16)	Head/Neck 16
Body coil (Body Matrix)	Body Matrix
Coil for peripheral angiography (PA Matrix).	PA Matrix
Knee coil (Extremity Matrix).	Extremity Matrix
Coil for lower body part (14-ch Extremity Coil).	14-ch Extremity Coil
Coil for foot (8 Ch Foot Ankle).	8 Ch Foot Ankle
Coil for wrist (8 Ch Wrist).	8 Ch Wrist
Coils for shoulders (Shoulder Array Small/Shoulder Array Large)	Shoulder Array Small/Large
Coils for breasts (Breast Matrix).	Breast Matrix
Endorectal coil (Endorectal).	Endorectal coil
Interface for flex coils connections	Connector
Flex large coils (Flex Large)	Flex Large
Flex small coils (Flex Small)	Flex Small
Coils for small parts (4 Ch Special Purpose)	Special purpose coil 4 Ch or 4 Ch Special Purpose

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