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Version 3 from 29.01.2020

High Fresh + solution, as part of:

- 1.1. High Fresh + solution in a bottle of 500 ml, 360 ml, 250 ml, 120 ml.
- 1.2. Lenses case (if necessary).

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Phone: (495) 937-2148

PRESCRIPTION

Multi-purpose solution “High Fresh+” is used to store, rinse, disinfect, remove protein and moisten all types of soft contact lenses.

INDICATIONS

Multi-purpose solution “High Fresh+” is recommended for rinsing and cleaning, disinfection and moistening soft contact lenses.

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DESCRIPTION

Multi-purpose solution High Fresh+ has an increased antibacterial and anti-protein effect. The pH of the solution is compatible with the eye tissue and organic eye fluids with which it contacts.

INSTRUCTIONS FOR USE

1. Wash your hands thoroughly.
2. Remove the lenses, gently wipe the lenses with a High Fresh+ solution on the palm of your hand for a few seconds.
3. Rinse the lenses with a High Fresh+ solution.
4. Place the lenses in a High Fresh+ solution for at least 4 hours.
5. Rinse the lenses with High Fresh+ solution before use.

WARNINGS AND PRECAUTIONS

After opening the package, use a High Fresh + solution for 60 days. Discard remaining unused solution.

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Multi-purpose solution High Fresh+ are supplied in bottles of various volumes: 500ml, 360ml, 250ml, 120ml, a pack with lens case or without lens case

CONTACT WITH HUMAN BODY

Soft contact lens care multi-purpose solutions "High Fresh+" have short-term contact with intact eye cornea of human. Solutions are washed off the lens with tear fluid and evaporated.

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USAGE CONDITIONS

Multi-purpose solution High Fresh+ is used at a temperature from +10 to +30 °C.

STORAGE CONDITIONS

Store at temperature from +5°C to +30°C

Relative air humidity from 35% to 70 %

Do not use after expiry date.

TRANSPORTATION CONDITIONS

Medical device can be transported by all means of transport. Multi-purpose solution "High Fresh+" is resistant to climate factors if transported at temperature from +5°C to +30°C and relative air humidity from 35% to 70 %.

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Medical device multi-purpose solution "High Fresh+" does not contain any medicinal drugs.

STERILIZATION

Multi-purpose solution "High Fresh+" in a bottle: 500ml, 360ml, 250ml, 120ml, a pack with lens case or without lens case is sterile.

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The bottle, tip of the bottle, cap, lens case for High Fresh+ solution, lens case, carton – not sterile.

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UTILIZATION

Multi-purpose solutions "High Fresh+", which have expired are to be disposed of in class A, as epidemiologically safe waste.

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Version 3 from 29.01.2020

High Fresh+ HYAL solution, as part of:

- 1.1. High Fresh + HYAL solution in a bottle of 360 ml, 250 ml, 120 ml, 20 ml.
- 1.2. Lenses case (if necessary).

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Phone: (495) 937-2148

PRESCRIPTION

Multi-purpose solution High Fresh+ HYAL is used to store, rinse, disinfect, remove protein and moisten all types of soft contact lenses.

INDICATIONS

Multi-purpose solution High Fresh+ HYAL is recommended for rinsing and cleaning, disinfection and moistening soft contact lenses.

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
DESCRIPTION

Multi-purpose solution High Fresh+ HYAL with sodium hyaluronate has an increased antibacterial and anti-protein effect. Due to the content of sodium hyaluronate, High Fresh+ HYAL solution is a unique balanced product that provides excellent comfort and long-term hydration of contact lenses.

INSTRUCTIONS FOR USE

1. Wash your hands thoroughly.
2. Remove the lenses, gently wipe the lenses with a High Fresh+ HYAL solution on the palm of your hand for a few seconds.
3. Rinse the lenses with a High Fresh+ HYAL solution.
4. Place the lenses in a High Fresh+ HYAL solution for at least 4 hours.
5. Before use, rinse the lenses with High Fresh+ HYAL solution.

WARNINGS AND PRECAUTIONS

After opening the package, use a High Fresh + HYAL solution for 60 days. Discard remaining unused solution.

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The multi-purpose High Fresh + HYAL solution is supplied in bottles of 360ml, 250ml, 120ml, 20ml, complete with a lenses case or without a lenses case.

CONTACT WITH HUMAN BODY

Soft contact lens care multi-purpose solutions High Fresh+ HYAL have short-term contact with intact eye cornea of human. Solutions are washed off the lens with tear fluid and evaporated.

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USAGE CONDITIONS

Multi-purpose solution High Fresh+ HYAL is used at a temperature from +10 to +30 °C.

STORAGE CONDITIONS

Store at temperature from +5°C to +30°C

Relative air humidity from 35% to 70 %

Do not use after expiry date.

TRANSPORTATION CONDITIONS

Medical device can be transported by all means of transport. Multi-purpose solution High Fresh+ HYAL is resistant to climate factors if transported at temperature from +5°C to +30°C and relative air humidity from 35% to 70 %.

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Medical device multi-purpose solution High Fresh+ HYAL does not contain any medicinal drugs.

STERILIZATION

Multi-purpose solution High Fresh+ HYAL in a bottle: 360ml, 250ml, 120ml, 20ml, a pack with lens case or without lens case is sterile.

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The bottle, tip of the bottle, cap, lens case for High Fresh HYAL+ solution, lens case, carton – not sterile.

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UTILIZATION

Multi-purpose solutions High Fresh+ HYAL, which have expired are to be disposed of in class A, as epidemiologically safe waste.

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A handwritten signature in black ink is written over a circular notary seal. The seal features a five-pointed star in the center and the text "LABATE FEDERICO" at the top and "PONTECORVO" at the bottom.

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Version 3 from 29.01.2020

High Fresh+ BIO solution, as part of:

1.1. High Fresh+ BIO solution in a bottle of 385 ml, 250 ml, 120 ml, 20 ml.

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Phone: (495) 937-2148

PRESCRIPTION

Multi-purpose solution High Fresh+ BIO is used to store, rinse, disinfect, remove protein and moisten all types of soft contact lenses.

INDICATIONS

Multi-purpose solution High Fresh+ BIO is recommended for rinsing and cleaning, disinfection and moistening soft contact lenses.

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DESCRIPTION

Multi-purpose solution High Fresh+ BIO has an increased antibacterial and anti-protein action. The solution contains sodium hyaluronate and plant-derived sources – extracts of Echinacea and sunflower seeds, which make its composition as natural as possible. Using a special extraction procedure, extracts from meristem of Echinacea stem cells are obtained. Echinacea stem cells are known for their moisturizing qualities, they contain soothing and moisturizing phytosterols and nourishing amino acids and minerals that provide osmotic balance. Sunflower lecithin extract, GPI Choline has soothing, moisturizing and nourishing effects due to the presence of nutrients such as choline, inositol and glycerin. Sodium hyaluronate polymer, with high viscoelastic qualities, improves lubrication and biocompatibility of lenses, maintaining their constant moisture.

INSTRUCTIONS FOR USE

1. Wash your hands thoroughly.

2. Remove the lenses, gently wipe the lenses with a High Fresh+ BIO solution on the palm of your hand for a few seconds.
3. Rinse the lenses with a High Fresh+BIO solution.
4. Place the lenses in a High Fresh+ BIO solution for at least 4 hours.
5. Before use, rinse the lenses with High Fresh+ BIO solution.

WARNINGS AND PRECAUTIONS

After opening the package, use a High Fresh + BIO solution for 60 days. Discard remaining unused solution.

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Do not remove the lid of the container while the disinfection process is in progress. After the disinfection process is completed, rinse the lens container and leave it to dry in the air.

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Multi-purposesolution High Fresh+ BIO in a bottle: 385ml, 250ml, 120ml, 20ml, a pack with lens case or without lens case

CONTACT WITH HUMAN BODY

Soft contact lens care multi-purpose solutions High Fresh+ BIO have short-term contact with intact eye cornea of human. Solutions are washed off the lens with tear fluid and evaporated.

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USAGE CONDITIONS

Multi-purpose solution High Fresh+ BIO is used at a temperature from +10 to +30 °C.

STORAGE CONDITIONS

Store at temperature from +5°C to +30°C

Relative air humidity from 35% to 70 %

Do not use after expiry date.

TRANSPORTATION CONDITIONS

Medical device can be transported by all means of transport. Multi-purpose solution High Fresh+ BIO is resistant to climate factors if transported at temperature from +5°C to +30°C and relative air humidity from 35% to 70 %.

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Medical device multi-purpose solution High Fresh+ BIO does not contain any medicinal drugs.

STERILIZATION

Multi-purpose solution High Fresh+ BIO in a bottle: 385ml, 250ml, 120ml, 20ml, a pack with lens case or without lens case is sterile.

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The bottle, tip of the bottle, cap, lens case for High Fresh+ BIO solution, lens case, carton – not sterile.

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UTILIZATION

Multi-purpose solutions High Fresh+ BIO, which have expired are to be disposed of in class A, as epidemiologically safe waste.

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WARRANTY POLICY

Company " Optoflexs.r.l." Italy, undertakes to guarantee to the end user that the medical device is free from defects and meets the manufacturer's specifications within 12 months from the date of sale, subject to the conditions specified in the instructions for use.

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INSTRUCTIONS FOR USE

Version 3 from 29.01.2020

High Fresh+ PEROXIDE solution, as part of:

- 1.1. High Fresh+ PEROXIDE solution in a bottle of 360 ml, 100 ml.
- 1.2. Lenses case.

Developer:

Optoflexs.r.l, Italy

Address: Via Aniense 10 – 00198 Rome (RM) Italy

Telephone: +39 06 4395508, Fax: + 39 06 4395486

E-mail: optoflex@legalmail.it

Place of manufacture:

Optoflexs.r.l, Italy

Address: Via Aniense 10 – 00198 Rome (RM) Italy

Authorized representative of the manufacturer:

“SVETLENZ” LLC

Address: 125363, Moscow, Fabritsiusast., build. 33, block 4, apt. 65

Phone: (495) 937-2148

PRESCRIPTION

Multi-purpose solution “High Fresh+” PEROXIDE is used to store, rinse, disinfect, remove protein and moisten all types of soft contact lenses.

INDICATIONS

Multi-purpose solution “High Fresh+” PEROXIDE is recommended for rinsing and cleaning, disinfection and moistening soft contact lenses.

CONTRAINDICATIONS

Do not use this product if you are allergic to any ingredient.

SIDE EFFECTS


If irritation or any other side effect occurs, remove the lenses and contact your eye care practitioner. Carefully follow the recommendations stated in the instructions for use.

DESCRIPTION

Multi-purpose solution High Fresh+ has an increased antibacterial and anti-protein effect. If used regularly, ensures complete sterilization of soft lenses.

INSTRUCTIONS FOR USE

1. Wash your hands thoroughly
2. Fill the case with multi-purpose solution High Fresh+ PEROXIDE up to the dose line. Use only the case with neutralizing platinum star that comes with the package of multi-purpose solution High Fresh+ PEROXIDE.
3. Lenses are placed in the lens storage case and immersed in the solution. The case must be kept upright. Providing that the solution functions properly, after the platinum star comes into the contact with the solution, oxygen bubbles will be generated. Otherwise, do not use the solution.

4. Store the lenses in the High Fresh+ PEROXIDE solution for at least 6 hours. During this period, when the hydrogen peroxide solution comes into reaction with the catalytic platinum star of the container, the lenses will be cleaned and disinfected, and the hydrogen peroxide will be completely neutralized.

5. Lenses can be worn after cleaning and disinfection. For particularly sensitive eyes, the lenses can be washed in saline before using.

WARNINGS AND PRECAUTIONS

After opening the package, use a high Fresh + solution for 60 days. Discard remaining unused solution.

Store solution at a temperature of +5 ° to +30 ° C.

Do not use solution after the expiration date.

Do not use solution if the cap is damaged before first use.

The solution is sterile while the cap is closed.

To avoid contaminating the solution do not touch dropper tip of bottle.

After use, firmly close the bottle with solution.

Do not remove the lid of the container while the disinfection process is in progress.

After the disinfection process is completed, rinse the lens container and leave it to dry in the air.

Always dispose of the solution left in the container after removing the lenses.

If the lenses have been in solution for more than 7 days, repeat the cleaning process.

Keep the solution out of the reach of children.

Never rinse contact lenses with tap water.

BASIC TECHNICAL CHARACTERISTICS

Name of material	Quantity	Manufacture
High Fresh+ PEROXIDE solution		
Hydrogen peroxide	3 %	Clean Consult, Italy
Phosphate buffer solution, stabilized ethylenediaminetetra acetic acid	0,03%	Clean Consult, Italy
Phosphonic acid	Osmolarity 290 - 320 mOsm/l	Clean Consult, Italy
pH	6,20 -6,90	

Multi-purpose solution "High Fresh+" PEROXIDE in a bottle: 360ml, 100ml, a pack with a lens case

CONTACT WITH HUMAN BODY

Soft contact lens care multi-purpose solutions "High Fresh+" PEROXIDE have short-term contact with intact eye cornea of human. Solutions are washed off the lens with tear fluid and evaporated.

Bottle, bottle tip, cap, lens case, carton have short-term (less than 24 hours) contact with intact skin.

USAGE CONDITIONS

Multi-purpose solutions High Fresh+ PEROXIDE are used at a temperature from +10 to +30 °C.

STORAGE CONDITIONS

Store at temperature from +5°C to +30°C

Relative air humidity from 35% to 70 %

Do not use after expiry date.

TRANSPORTATION CONDITIONS

Medical device can be transported by all means of transport. Multi-purpose solution "High Fresh+" PEROXIDE is resistant to climate factors if transported at temperature from +5°C to +30°C and relative air humidity from 35% to 70 %.

INFORMATION ON COMPLIANCE WITH THE DATA OF THE STATE REGISTER OF MEDICINES FOR MEDICAL USE CONTAINED IN THE MEDICAL PRODUCT

Medical device multi-purpose solution "High Fresh+" PEROXIDE does not contain any medicinal drugs.

STERILIZATION

Multi-purpose solution "High Fresh+" PEROXIDE in a bottle: 360ml, 100ml, a pack with lens case.

Sterilization using aseptic treatment methods.


The bottle, tip of the bottle, cap, lens case for High Fresh+ PEROXIDE solution, lens case, carton – not sterile..

SHELF LIFE

Shelf life of solution in undamaged packaging is 3 years from the date of manufacture. Useful life of solution is 60 (sixty) days after opening the bottle.

UTILIZATION

Multi-purpose solutions "High Fresh+" PEROXIDE, which have expired are to be disposed of in class A, as epidemiologically safe waste.

After using the container (bottle), as well as the case for lenses are to be disposed of in class A, as epidemiologically safe waste. Therefore, disposal of used products must be carried out in accordance with good medical practice, as well as applicable laws and regulations. Disposal must be carried out in accordance with the rules and regulations approved by Russian law

LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS

Medical device conforms to the requirements of the following standards:

Directive 93/42/EEC «Medical tools, devices, equipment»

EN ISO 13485 «Medical devices. Quality management systems. Requirements for regulatory purposes»

ISO 10993-5 «Medical devices. Biological evaluation of medical devices. Part 5. Tests for *in vitro* cytotoxicity»

EN ISO 14729 – 2010 «Ophthalmic optics. Contact lens care products. Microbiological requirements and test methods for products. Regimens for hygienic management of contact lenses».

EN ISO 14730- 2000 «Ophthalmic optics. Contact lens care products. Antimicrobial preservative efficacy testing and guidance on determining discard date».

WARRANTY POLICY

Company "Optoflexs.r.l." Italy, undertakes to guarantee to the end user that the medical device is free from defects and meets the manufacturer's specifications within 12 months from the date of sale, subject to the conditions specified in the instructions for use.

If the product is found to be defective before the end of the warranty period, your exclusive right and sole obligation to "Optoflexs.r.l.", Italy is to replace faulty products of the brand "Optoflexs.r.l." Italy.

For any questions related to the use of the product, please contact customer support or contact your local distributor.

Manufacturer:

Optoflexs.r.l., Italy

Address: Via Aniene 10 – 00198 Rome (RM) Italy

Telephone: +39 06 4395508, Fax: + 39 06 4395486

E-mail: optoflex@legalmail.it

Regional distributor in Russia:

“SVETLENZ” LLC

Address: 125363, Moscow, Fabritsiusast., build. 33, block 4, apt. 65

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