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Apostille

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1. Country (Pays): **United States of America**

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4. bears the seal / stamp of (est revêtu du sceau / timbre de)
Ryan J Bartz, Notary Public, State of Minnesota

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
To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Operator Manual for 3M™ Attest EO Biological Indicators Incubator in English is truthful and accurate.

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
Operator's Manual

Attest Biological Indicator Incubators Attest EO Biological Indicators Incubator

Product name

Attest EO Biological Indicators Incubator.

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The 3M Attest EO Incubator is designed and has been tested to be used with the 3M Attest Biological Indicators (BI) 1264 and 1264P and allows rapid, reliable monitoring of ethylene oxide sterilization process. As a quick reference, the ethylene oxide monitoring system is color-coded green. Use with other biological indicators may lead to incorrect results.

Each Attest incubator is pre-set at the appropriate temperature to promote growth of the test organism. The incubator's internal temperatures are rigidly controlled to 37°C +/-2°C (99°F +/-3°F) and automatically maintained. If a temperature measurement device is used to monitor the incubator, the sensor should be in contact with the metal heating block wall.

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SAFETY INFORMATION

Warnings

- To reduce the risk of electric shock:
 - do not pour any liquids into the metal heating block holes
 - do not immerse the incubator in any liquid
 - use indoors only
 - disconnect power cord before cleaning
 - use only a disinfectant or mild detergent and wipe with a damp cloth
 - there are no user serviceable parts
 - the incubator must be returned to the manufacturer for repair
 - make sure a properly rated power cord is installed on the incubator.
- To reduce the risk of using a load from a failed sterilization cycle:
 - match the color of the BI cap with the same color as the stripe on the incubator metal heater block
 - verify the incubation conditions are correct by using a positive control BI each day a test BI is incubated.
- There is a glass ampoule inside the plastic vial of the biological indicator:
 - crushing or excessive handling of the BI before cooling may cause the glass ampoule to burst
 - wear safety glasses and gloves when removing the BI from the sterilizer
 - wear safety glasses when crushing the biological indicator
 - handle the BI by the cap when crushing
 - do not use your fingers to crush the glass ampoule.

Precautions

- To avoid possible environmental contamination, dispose of the incubator according to applicable governmental regulations.
- Avoid contact with the hot surface of the metal heating block.
- If this product is used in a manner not specified, the protection provided by the product may be impaired.

Explanation of Product and Packaging Label Symbols (Note: not all symbols appear on all models)

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Manufacturer

Authorized Representative for the European Community

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Contents of delivery

Biological Indicator Incubator with cap and power cord – 1 PC.

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Environmental Condition	Operating Condition
Altitude	2500 (max) meters

Operating Temperature	15 – 35 °C
Relative Humidity	35 – 80% (non-condensing)
Voltage Range	220-240 Volts AC
Frequency	50/60 Hertz
Current	0.2 Amps
Degrees of protection provided by enclosures (IP Code)	IPX0 (no moisture protection inside the housing)
Electrical safety class	Class I
Installation/Over Voltage	Category II
Pollution Degree	2

SPECIFICATIONS*

* The tolerance is $\pm 5\%$ for all values, unless otherwise specified

Parameter/characteristic	129
Dimensions, mm	275 x 170 x 125 mm
Weight, kg	No more than 2
The number of cells for incubating indicators	28
Incubation temperature	37°C $\pm$ 2
Incubation time	12-24 h (preliminary result) 48 h (final result)
Application area	Ethylene oxide sterilization
Attest Biological Indicators	1264 & 1264P
Color Coding	Green
The function of self-diagnosis of the internal temperature of the incubator	+
Function of automatic reading of biological indicators (negative fluorescence)	–
Power cord length, cm	230

Directions for Use

1. Plug incubator in and allow to warm-up before use. To eliminate warm-up periods, it is recommended that the unit be left on and that the cover remains closed except during BI insertion and removal. The LED lights in the housing indicate when the incubator heater blocks have reached the correct temperature. The LED indicator on the left corresponds to the heater block in the upper tier and the LED indicator on the right corresponds to the heater block in the lower tier. A red light indicates the unit is powered on but the heater block temperature is below the set point. Do not incubate biological indicators until the LED indicators turn green, indicating the incubator has reached the temperature set point.
2. Remove the biological indicator test pack from load according to facility policy. Remove the Attest biological indicator from the test pack according to the product instructions for use.
3. Position the biological indicator in the incubator's metal heating block (See Figure 1). Place the bottom of the indicator into the metal heating block so that the indicator is at an angle of approximately 45° and the cap is above the metal block.
4. Push the indicator straight back (See Figure 2). This crushes the ampoule and activates the indicator.
5. Push the activated indicator down to seat it in the metal heating block (See Figure 3). Be sure that the cap remains above the metal block. Refer to the Attest Biological Indicator 1264 Instructions for Use for further directions for use, including the interpretation of results and use of controls.

Connectors

- IEC C14 connector for connecting the power cord supplied by «3M».


Figure 1


Figure 2


Figure 3

- **Storage and Transportation:** Temperature from -10°C to +50°C, Relative Humidity Limitation from 20% to 80% (non-condensing), a pressure are not regulated.
During transportation the products may experience short term temperature variations without affecting product safety and performance.
Transportation of medical device is carried out by all type of transport.

Cleaning and disinfection

Clean the exterior surface of the 3M Attest Biological Indicator Incubators by wiping it with a damp cloth moistened with a mild detergent. Do not attempt to clean interior components. If additional cleaning of the wells is required, call for repair or replacement.

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Electromagnetic Compatibility – 220V Incubators – 28 Vial

The Biological Incubators completed conducted and radiated emissions testing and immunity testing in accordance with Directive 2004/108/EC - Electromagnetic Compatibility Directive (EMC). These models were evaluated to the following European Standard EN 61326-1 (2006); IEC 61326-1 (2005) requirements listed below and found to be compatible with its electromagnetic environment.

EN/IEC 61326-1 Radiated/Conducted Emissions

EN/IEC 61000-4-3 Radiated Immunity

EN/IEC 61000-3-2 Harmonic Emissions

EN/IEC 61000-3-3 Flicker Emissions

EN/IEC 61000-4-2 Electrostatic Discharge (ESD)

EN/IEC 61000-4-4 Electrical Fast Transient (EFT)

EN/IEC 61000-4-5 Surge

EN/IEC 61000-4-6 Conducted Immunity

EN/IEC 61000-4-11 Voltage Dips and Interruption

EN/IEC 61000-4-8 Power Frequency Magnetic Fields

Disposal

Dispose of content/container in accordance with the local/regional/national/international/regulations.

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The product, as well as other accessories, must not be disposed of with unsorted household waste.

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Warranty

The 3M Attest EO Incubator has a one year limited product warranty. The US warranty, remedy and limitations are described in the warranty card provided with the device. For countries outside of the US, any warranty is established by the subsidiary with responsibility for service of the incubator.

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Attest Biological Indicator Incubators
Attest Steam Biological Indicators Incubator

Product name

Attest Steam Biological Indicators Incubator.

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The 3M Attest Steam Incubator is designed and has been tested to be used with 3M Attest Biological Indicators (BI) 1261 and 1261P or 1262 and 1262P and allows rapid, reliable monitoring of steam sterilization processes. As a quick reference, the steam monitoring system is color-coded blue and brown, respectively. Use with other biological indicators may lead to incorrect results.

Each Attest incubator is pre-set at the appropriate temperature to promote growth of the test organism. The incubator's internal temperatures are rigidly controlled to $56^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($133^{\circ}\text{F} \pm 3^{\circ}\text{F}$) and automatically maintained. If a temperature measurement device is used to monitor the incubator, the sensor should be in contact with the metal heating block wall.

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SAFETY INFORMATION

Warnings

- To reduce the risk of electric shock:
 - do not pour any liquids into the metal heating block holes
 - do not immerse the incubator in any liquid
 - use indoors only
 - disconnect power cord before cleaning
 - use only a disinfectant or mild detergent and wipe with a damp cloth
 - there are no user serviceable parts
 - the incubator must be returned to the manufacturer for repair
 - make sure a properly rated power cord is installed on the incubator.
- To reduce the risk of using a load from a failed sterilization cycle:
 - match the color on the exterior label of the incubator with the same color cap on the BI
 - verify the incubation conditions are correct by using a positive control BI each day a test BI is incubated.
- There is a glass ampoule inside the plastic vial of the biological indicator:
 - crushing or excessive handling of the BI before cooling may cause the glass ampoule to burst
 - wear safety glasses and gloves when removing the BI from the sterilizer
 - wear safety glasses when crushing the biological indicator
 - handle the BI by the cap when crushing
 - do not use your fingers to crush the glass ampoule.

Precautions

- To avoid possible environmental contamination, dispose of the incubator according to applicable governmental regulations.
- Avoid contact with the hot surface of the metal heating block.
- If this product is used in a manner not specified, the protection provided by the product may be impaired.

Explanation of Product and Packaging Label Symbols (Note: not all symbols appear on all models)

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Manufacturer

Authorized Representative for the European Community

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Contents of delivery

Biological Indicator Incubator with cap and power cord - 1 PC.
Operator's manual - 1 PC.

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Environmental Condition	Operating Condition
Altitude	2500 (max) meters
Operating Temperature	15 - 35 °C
Relative Humidity	20 - 80% (non-condensing)
Voltage Range	220-240 Volts AC

Frequency	50/60 Hertz
Current	0.13 Amps
Degrees of protection provided by enclosures (IP Code)	IPX0 (no moisture protection inside the housing)
Electrical safety class	Class I
Installation/Over Voltage	Category II
Pollution Degree	2

Attest Incubator 220-240 volt	Sterilization Process Monitored	Attest Biological Indicators	Color Coding
14 vials 56 +/- 2°C (133 +/- 3°F)	Steam (except 270°F/132°C gravity displacement)	1262 & 1262P	Brown
	Steam gravity displacement 270°F/132°C	1261 & 1261P	Blue

SPECIFICATIONS*

* The tolerance is $\pm 5\%$ for all values, unless otherwise specified


Parameter/characteristic	118
Dimensions, mm	100 x 100 x 130
Weight, kg	No more than 1
The number of cells for incubating indicators	14
Incubation temperature	56°C $\pm$ 2
Incubation time	12-24 h (preliminary result) 48 h (final result)
Application area	Steam sterilization
The function of self-diagnosis of the internal temperature of the incubator	+
Function of automatic reading of biological indicators (negative fluorescence)	-
Power cord length, cm	230

Directions for Use

1. Plug incubator in and allow to warm-up for 30 minutes before use. It is recommended that the unit be left on to eliminate warm-up periods.
2. After completion of the cycle and while wearing safety glasses and gloves, fully open the sterilizer door for a minimum of 5 minutes prior to removing the Attest biological indicator.
3. When the biological indicator is not contained in a test pack or other heat absorbing packaging material, remove the biological indicator from the sterilizer and allow to cool for an additional 10 minutes prior to crushing.
4. When the biological indicator is contained in a test pack or other heat absorbing packaging material, the test pack or any other heat absorbing packaging material should be removed from the sterilizer and opened up for 5 minutes to dissipate heat prior to removing the biological indicator. Then allow the indicator to cool outside the test pack for an additional 10 minutes prior to crushing.
5. Position the biological indicator in the incubator's metal heating block (See Figure 1). Place the bottom of the indicator into the metal heating block so that the indicator is at an angle of approximately 45° and the cap is above the metal block.
6. Push the indicator straight back (See Figure 2). This crushes the ampoule and activates the indicator.
7. Push the activated indicator down to seat it in the metal heating block (See Figure 3). Be sure that the cap remains above the metal block.

Refer to the Attest biological indicator Instructions for Use for further directions for use, including the interpretation of results and use of controls.

Storage and Transportation: Temperature from -10°C to +50°C, Relative Humidity

to 80% (non-condensing), a pressure are not regulated.

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Cleaning and disinfection

Clean the exterior surface of the 3M Attest Biological Indicator Incubators by wiping it with a damp cloth moistened with a mild detergent. Do not attempt to clean interior components. If additional cleaning of the wells is required, call for repair or replacement.

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Electromagnetic Compatibility – 220V Incubators – 14 Vial

The Biological Incubators completed conducted and radiated emissions testing and immunity testing in accordance with Directive 2004/108/EC - Electromagnetic Compatibility Directive (EMC). These models were evaluated to the following European Standard EN 61326-1 (2006); IEC 61326-1 (2005) requirements listed below and found to be compatible with its electromagnetic environment.

EN/IEC 61326-1 Radiated/Conducted Emissions

EN/IEC 61000-4-3 Radiated Immunity

EN/IEC 61000-3-2 Harmonic Emissions

EN/IEC 61000-3-3 Flicker Emissions

EN/IEC 61000-4-2 Electrostatic Discharge (ESD)

EN/IEC 61000-4-4 Electrical Fast Transient (EFT)

EN/IEC 61000-4-5 Surge

EN/IEC 61000-4-6 Conducted Immunity

EN/IEC 61000-4-11 Voltage Dips and Interruption

EN/IEC 61000-4-8 Power Frequency Magnetic Fields

Disposal

Dispose of content/container in accordance with the local/regional/national/international/regulations.

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- The WEEE logo on the packaging, the appliance and in the instructions indicates that disposal must be carried out in accordance with the EU Electronic and Electrical Equipment Directive.

The product, as well as other accessories, must not be disposed of with unsorted household waste.

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AO «3M Russia»

108811, Moscow, p. Moskovsky, Kievskoe highway, 22nd km, house. 6, b. 1
Telephone numbers: +7 (495) 784 7474 / +7 (495) 784 7475 (fax)
e-mail: 3mra@mmm.com

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
To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Operator Manual for 3M™ Attest Rapid Readout Steam Biological Indicators Incubator (Auto-reader)

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

Attest Biological Indicator Incubators
Attest Rapid Readout Steam Biological Indicators Incubator
(Auto-reader)

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Hardware Disclaimer

The 3M Attest Auto-reader 390 hardware specifications are subject to change. The system images, hardware components and hardware specifications included in the manual may not match the system installed. Any changes or modifications to the authorized system installation have been verified as compatible with the functionality outlined in this document.

Product name

Attest Rapid Readout Steam Biological Indicators Incubator (Auto-reader)

(Hereinafter, the devices may be referred to as "Biological Indicators Incubator (etc.), "Incubator" and etc.))

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The 3M Attest Auto-reader 390 is designed to incubate and automatically read the 3M Attest Rapid Readout Biological Indicators for Steam, 1291 and 1292, at 60°C for a final fluorescent result at 1 hour for 1291 and 3 hours for 1292.

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No known possible side effects.

Product description

The 3M Attest Rapid Readout Biological Indicator System provides rapid, reliable monitoring of steam sterilization processes performed in health care facilities. The system is comprised of biological indicators, challenge packs, and an incubator/reader (3M Attest Auto-reader 390) with embedded software.

The 3M Attest Auto-reader 390 is designed to incubate and automatically read 3M Attest 1291 Rapid Readout Biological Indicators for Steam at 60°C for a final fluorescence result at 1 hour and 3M Attest 1292 and 1292E Rapid Readout Biological Indicators for Steam at 60°C for a final fluorescence result at 3 hours. The 3M Attest 1292E Rapid Readout Biological Indicator is not cleared for use in the United States. The 3M Attest Auto-reader 390 is also designed to allow for further incubation of 3M Attest Rapid Readout Biological Indicators for Steam, catalog numbers 1291 and 1292, for a final negative visual pH color change result at 24 hours and 48 hours, respectively. This further incubation is optional and not intended for routine use.

A positive fluorescence reading or a visual pH color change indicates a steam sterilization process failure. However, due to the high sensitivity of the fluorescent result, there is no advantage to continuing to incubate processed 3M Attest Rapid Readout Biological Indicators once the fluorescent result has been documented. Other biological indicators (BIs) are not compatible with this device and cannot be used.

The unit is designed to be used only with the Power Supply module and Ethernet cable supplied by 3M.

Serial Number

For easy identification, each 3M Attest Auto-reader 390 has a unique serial number printed on a label found on the back of the unit and displayed on the left side of the upper row of text on the LCD Display when the user depresses the  button. Please record your serial number in this manual for future reference:

SAFETY AND PRECAUTIONS

The 3M Attest Auto-reader 390 and its related devices and accessories are designed to provide safe and reliable service when used according to the instructions provided. Please read, understand, and follow all safety information contained in the instructions for use included with the 3M Attest Auto-reader 390 and Rapid Readout Biological Indicator devices prior to use. Use this equipment only for the purpose described in this Operator's Manual. **Retain these instructions for future reference.**

If this product is used in a manner not specified, the protection provided by the product may be impaired.

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Always unplug the 3M Attest Auto-reader 390 and allow to cool before cleaning.

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Do not open the instrument housing – there are no user serviceable parts.

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Allow the 3M Attest Rapid Readout Biological Indicator to cool for the recommended time period before activating. Activating or excessive handling of the biological indicator before cooling may cause the glass ampoule to burst.

Wear safety glasses when activating the 3M Attest Rapid Readout Biological Indicator.

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Instrument is to be used by operators familiar with the device, device functionality, and Operator's Manual.

Do not place the instrument in environment exposed to sunlight or strong incandescent light.

Do not place instrument close to any device that emits a strong electro-magnetic field.

Do not use on a vibrating surface.

Match the cap color of the 3M Attest Rapid Readout Biological Indicator with the color-coded configuration sticker surrounding the 3M Attest Auto-reader 390 incubation wells.

Do not remove or change placement of 3M Attest Rapid Readout Biological Indicator once it is placed

into a well.

Do not remove the 3M Attest Rapid Readout Biological Indicator from the incubation well until the (+) or (-) symbol on the LCD panel indicates the test is complete.

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HEALTH AND SAFETY INFORMATION

Device Safety Compliance

The 3M Attest Auto-reader 390 complies with the following standards as demonstrated by the CB Scheme Certificate and test report issued by Underwriters Laboratories (UL):

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The 3M Attest Auto-reader 390 is listed as Laboratory equipment and carries the UL mark with adjacent indicators "C" and "US" based on compliance to the standards UL 61010-1 and CAN/CSA 22.2 No. 61010-1.

The 3M Attest Auto-reader 390 complies with the CE mark related to the Low Voltage Directive (LVD) 2006/95/EC as confirmed in the Declaration of Conformity.

The 3M Attest Auto-reader 390 complies with the RoHS Directive, Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

The 3M Attest Auto-reader 390 complies with the WEEE Directive, Directive 2012/19/EU of the European Parliament and of the Council of 04 July 2012 on waste electrical and electronic equipment (WEEE).

EMC Compliance

The 3M Attest Auto-reader 390 complies with the following EMC standards as confirmed in the Certificate of Compliance generated by 3M:

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The 3M Attest Auto-reader 390 complies with the EMC requirements of the CE mark EMC Directive 2004/108/EC.

The 3M Attest Auto-reader 390 complies with the Australian and New Zealand electrical safety and electromagnetic compatibility requirements as confirmed in the Supplier's Declaration of Conformity that is linked to the Australian/New Zealand RCM (Regulatory Compliance Mark).

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Class 1 LED Product

An investigation was conducted on the 3M Attest Auto-reader 390 per the requirements of IEC 60825-1 Ed. 1.2 (A2:2001). The 3M Attest Auto-reader 390 is a Class 1 LED product with Class 1 Internal Radiation.

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Note: The use of the term "controls" in the above caution statement as required by IEC 60825-1 Ed. 1.2 (A2:2001) refers to mechanisms by which the 3M Attest Auto-reader 390 may be controlled and not to the use of biological indicator controls.

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
The 3M Attest Auto-reader 390 is not designed to incubate 3M Attest 1261 and 1262 Steam Biological Indicators, which have a steam entry port in the center of the top of the cap. The 3M Attest Auto-reader 390 is designed to incubate and read 3M Attest Rapid Readout Biological Indicators for Steam, catalog numbers 1291, 1292, and 1292E.

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Contents of delivery

1. Biological Indicator Incubator – 1 PC.
2. Color-coded configuration sticker – 3 PC.
3. Power supply – 1 PC.
4. Power adapter – 4 PC.
5. Ethernet cable – 1 PC.
6. Operator's manual – 1 PC.

SPECIFICATIONS*

* The tolerance is $\pm 5\%$ for all values, unless otherwise specified

Parameter/characteristic	390
Dimensions, mm	350 x 180 x 130
Weight, kg	No more than 2
The number of cells for incubating indicators	10
Incubation temperature	$60^{\circ}\text{C} \pm 2^{\circ}\text{C}$
Incubation time	1 h (for biological indicator 132°C , 1291) 3 hours (for biological indicator 121°C , 1292) If additional research is needed: incubation - 7 hours, determination of the color change of the pH indicator
Application area	Steam sterilization
The function of self-diagnosis of the internal temperature of the incubator	+
Function of automatic reading of biological indicators (negative fluorescence)	+
Power Supply Specifications	
Input characteristics (power supply)	

Parameter/characteristic	390
Voltage Range	100-240 Volts AC
Frequency	50/60 Hertz
Current	0,6 Amps
Output characteristics (power supply)	
Voltage Range	12 Volts DC
Current	1,5 Amps
Power supply cord length	150 cm
Battery	
Type	CR 2032
Max reverse charging current	10mA
Software	
Software name	390
Carrier	3M Company
Functionality	Incubation of BI and measurement
Version	3.0.6
Release date	June 16, 2015
Classification (in accordance with IEC 62304)	Class B
Volume	The code size of the master PIC 128k bytes and the size of the RAM 3808 bytes

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Environmental Condition	Operating Condition
Altitude	3000 (max) Meters
Operating Temperature	16 – 40 °C
Operating Relative Humidity	20 – 80% (non-condensing)
Voltage Range	12 Volts DC
Frequency	DC
Current	1,2 Amps
Degrees of protection provided by enclosures (IP Code)	IPX0 (no moisture protection inside the housing)
Electrical safety class	Class I
Installation/Over Voltage	Category II
Pollution Degree	2

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Choose Location

Place 3M Attest Auto-reader 390 on a firm level surface away from sunlight and incandescent light. To save counter space, units may be stacked two high.

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Select & Install Correct Power Plug

The 3M Attest Auto-reader 390 power supply comes with multiple blade configurations. Install the blade configuration appropriate for your country into the power supply and plug into appropriate AC grounded outlet.

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Wait 30 Minutes Before Using


Once plugged in, allow a 30 minute warm-up period before placing 3M Attest Rapid Readout Biological Indicators into the incubation wells. The C1 caution code will disappear when proper incubation

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INSTRUCTIONS FOR USE – Wells, Buttons, & Display Interpretation

Please familiarize yourself with the face of the 3M Attest Auto-reader 390:

Wells

1. Incubation/Reader Wells – The 3M Attest Auto-reader 390 has ten incubation wells.

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10. LCD Display – Alarm Enabled Indicator

Color-Coded Configuration Indicator

11. Color-coded incubation well configuration sticker

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INSTRUCTIONS FOR USE – 3M Attest Auto-reader 390 Configuration

Note: Ideally, the 3M Attest Auto-reader 390 should be empty (i.e., have no incubating BIs) when it is reconfigured.

Press the  and  buttons on the front of the 3M Attest Auto-reader 390 simultaneously for 3 seconds to enter configuration mode. Once in configuration mode, use the Up  and Down  buttons to set the

3M Attest Auto-reader 390 clock, audible alarm, and incubation well configuration as desired.
The button  acts as an "Enter" key and moves the flashing cursor to the next changeable setting.

Setting the 3M Attest Auto-reader 390 Clock

The 3M Attest Auto-reader 390 date and time are displayed in this format: YYYY-MM-DD HH:MM, with the time displayed using a 24 hour clock. To change the year, if necessary, push the up or down buttons to set desired year, then press  to advance to the month, adjust as necessary and press  to advance to the day. Use the  button to move to the time and set the hours and minutes.

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Note: The 3M Attest Auto-reader 390 display backlight flashes whenever an audible alarm sounds. Disabling the audible alarm does NOT disable the LCD display backlight flashing.

Press  to select the desired option and advance to the bottom row of text.

Setting the Incubation Well Configuration

Incubation wells 1-5 and 6-10 can be configured to run the same or different types of 3M Attest Rapid Readout Biological Indicators. The flashing "A" in Attest indicates the setting may be changed by the user. There are three incubation well configuration options:

LCD Panel Configuration Setting		Result
Attest 1291	Attest 1292	Auto-reader wells 1 - 5 configured for 1291 BIs (blue cap, 1 hour readout) and wells 6 - 10 configured for 1292 BIs (brown cap, 3 hour readout)
Attest 1291	Attest 1291	Auto-reader wells 1 -10 configured for 1291 BIs (blue cap, 1 hour readout)
Attest 1292	Attest 1292	Auto-reader wells 1 -10 configured for 1292 BIs (brown cap, 3 hour readout)

Use the Up Arrow or Down Arrow buttons to move through the three options until the desired configuration setting is displayed.

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Note 1: If there are any incubating BIs which do not yet have a final result (+ or -) displayed, the well configuration for the occupied wells will NOT change until the final result is displayed on the LCD panel and the BI is removed from the reader.

Note 2: Ensure the color-coded incubation well configuration sticker surrounding the incubation wells matches the programmed well configuration.

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
WARNING: Match the cap color of the 3M Attest Rapid Readout Biological Indicator with the color-coded configuration sticker surrounding the 3M Attest Auto-reader 390 incubation wells.

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 **CAUTION:** Allow the 3M Attest Rapid Readout Biological Indicator to cool for the recommended time period before activating. Activating or excessive handling of the biological indicator before cooling may cause the glass ampoule to burst.

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
Follow this procedure for every 3M Attest Rapid Readout Biological Indicator to be read.


While wearing safety glasses, **CLOSE** the cap of the cooled 3M Attest Rapid Readout Biological Indicator (BI) by pushing cap down.

CRUSH the glass media ampoule in the Crusher Well built into the top of the 390 Auto-reader.

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PLACE the crushed BI in **color-coded** incubation well. Observe 3M Attest Auto-reader 390 LCD panel to confirm correct placement of the BI. If the BI is correctly seated, the 3M Attest Auto-reader 390 LCD panel will display ***, **, * and then remaining minutes of incubation directly below the incubating BI.
RECORD BI result.

 **WARNING:** Do not remove the 3M Attest Rapid Readout Biological Indicator from the incubation well until the (+) or (-) symbol on the LCD panel indicates the test is complete.

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BI Placement

Wells must be configured and color coded to match the cap of the 3M Attest Rapid Readout Biological Indicator, e.g. the 3M Attest 1291 Rapid Readout Biological Indicator with a blue cap goes into a blue well and the 3M Attest 1292 Rapid Readout Biological Indicator with a brown cap goes into a brown well.

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Note: If 3M Attest Rapid Readout Biological Indicator is inadvertently removed before the incubation period is complete, an alarm will sound and a C2 Caution code will be displayed on the LCD panel under that well. You must return the 3M Attest Rapid Readout Biological Indicator to the reader well within 10 seconds or the results will be lost.

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Positive Control

Activate and incubate at least one non-processed 3M Attest Rapid Readout Biological Indicator to use as the positive control each day that a processed 3M Attest Rapid Readout Biological Indicator is incubated and read by the unit. Write a "C" (for "control") and the date on the label. The positive control can be placed in any appropriately configured (i.e., blue or brown) incubator well. The positive control should be the same part number and have the same lot number as the processed 3M Attest Rapid Readout Biological Indicator in the unit.

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- correct incubation temperatures are met (heater block reading of 60°C +/-2°C);
- viability of spores has not been altered due to improper storage temperature, humidity, or proximity to chemicals;
- capability of media to promote rapid growth; and
- proper functioning of 3M Attest Auto-reader 390 components.

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Fluorescent Results

After placing a 3M Attest Rapid Readout Biological Indicator into a corresponding color-coded well, the 3M Attest Auto-reader 390 LCD panel will display ***, **, * and then the remaining minutes of incubation directly below the incubating BI to indicate that incubation/reading is in progress.

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At the end of the specified incubation time, if a negative 3M Attest Rapid Readout Biological Indicator result is detected, a (-) will display on the LCD panel directly under the incubating BI. The final negative fluorescent result (-) for a processed 3M Attest Rapid Readout Biological Indicator after 1 hour or 3 hours of incubation, for 1291 and 1292 BIs, respectively, indicates an acceptable sterilization process. The positive control should read positive (+).

After either the + or - symbol is displayed, and the result is recorded, the 3M Attest Rapid Readout Biological Indicator can be removed from the well. The 3M Attest Rapid Readout Biological Indicator should not be removed or touched when the display indicates remaining incubation time.

The positive (unprocessed) control 3M Attest Rapid Readout Biological Indicator must provide a positive fluorescent result (+). Processed 3M Attest Rapid Readout Biological Indicator results are not valid until the positive control reads fluorescent positive (+). With a processed 3M Attest Rapid Readout Biological Indicator, a positive (+) result indicates a steam sterilization process failure.

Act on any positive test results immediately. Determine the cause of the positive 3M Attest Rapid Readout Biological Indicator following current facility policies and procedures. Always retest the sterilizer and do not use the sterilizer for processing loads until 3 consecutive 3M Attest Rapid Readout Biological Indicator results are negative.

Optional Visual pH Color Change Result

3M Attest Rapid Readout Biological Indicators may be further incubated for a visual pH color change. This optional further incubation is not intended for routine use and may be used for special studies. Due to the high sensitivity of the fluorescent result, there is no advantage to incubating the 3M Attest Rapid Readout Biological Indicators 1291 and 1292 beyond 1 and 3 hours, respectively. In the case of a positive control BI, a yellow color change of the growth media will appear within 24 hours for 1291 BIs and within 48 hours for 1292 BIs. Any observation of a yellow color within the vial indicates a positive result. In the case of a processed BI, a media color change from purple to yellow indicates a steam sterilization process failure. A negative pH color change result, i.e., media remains purple, can be assessed at 24 hours for 1291 BIs and at 48 hours for 1292 BIs.

FEATURES OF THE 3M ATTEST AUTO-READER 390

Audible Alarm

An audible alarm will sound, the LCD display backlight will flash, and a (+) symbol will be displayed on the LCD panel directly under the incubating BI when a positive result is detected. An alarm will also sound if a 3M Attest Rapid Readout Biological Indicator is removed from the 3M Attest Auto-reader 390 well before a final result has been detected. Data will be lost if the 3M Attest Rapid Readout Biological Indicator is not replaced in the incubation well within 10 seconds.

Alarm On/Off

The audible alarm for positive results can be turned off for a particular 3M Attest Rapid Readout Biological Indicator. This is accomplished by pressing the "Buzzer Alarm OFF" button after the alarm has sounded. This will turn off the alarm for that BI only. Any positive 3M Attest Rapid Readout Biological Indicator results subsequently detected in other wells will result in the alarm sounding. Removing the 3M Attest Rapid Readout Biological Indicator will reactivate the alarm for that well. Pressing the "Buzzer Alarm OFF" button will not silence the alarm for error or caution code alarm warnings. To disable the alarm, see configuration section.

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Stackable

To save counter space, 3M Attest Auto-reader 390s may be stacked two high.

Information Button

When the Information button is pressed once the 3M Attest Auto-reader 390 displays two lines of text in the LCD display.

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Press the information button twice and the 3M Attest Auto-reader 390 LCD panel will display two additional lines of text.

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- Date (YYYY-MM-DD) and time that power was supplied to the 3M Attest Auto-reader 390.

Temperature Monitoring

The 3M Attest Auto-reader 390 contains an internal temperature self-diagnostic program. Should the

temperature fall outside the specified range (60 °C +/-2°C) an error code will be displayed. The heater block temperature may be checked by depressing the Temperature button once.

Easy Configuration

The 3M Attest Auto-reader 390 may be configured using either the buttons located on the front of the unit or using the optional 3M Attest Auto-reader 390 Web App.

Optional 3M Attest Auto-reader 390 Web App

If desired, the 3M Attest Auto-reader 390 may be connected to a PC via an Ethernet connection for additional capabilities:

- Remotely view current status of BIs in the 3M Attest Auto-reader 390 via PC or mobile hand held devices. For full functionality, the use of an HTML5 compliant browser with JavaScript enabled is recommended.
- Remotely view and configure 3M Attest Auto-reader 390 settings.
- View the last 100 BI records.
- Print out documentation using a local or networked printer and/or save documentation as an electronic file.*

Note: To access the 3M Attest Auto-reader 390 Web App within a secure network, the PC or mobile device must have access to the secure network; options include connecting via a network Ethernet cable, signing in to the facility's wireless network, or connecting via a virtual private network (VPN).

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CARE AND CLEANING

	CAUTION: Always unplug the 3M Attest Auto-reader 390 and allow to cool before cleaning.
	Do not spill liquid onto or into the instrument. Do not immerse the unit in liquid.
	Clean external surfaces using only the instructions provided by the manufacturer.
	The instrument must be returned to the manufacturer for repair.

Clean the exterior surface of the 3M Attest Auto-reader 390 by wiping it with a damp cloth moistened with a mild detergent. Do not attempt to clean interior components. If additional cleaning of the wells is required, call for repair or replacement.

Note: The 3M Attest Auto-reader 390 does not contain any user serviceable parts or require any routine maintenance. Any malfunctions will be indicated by an error code. Refer to the Troubleshooting Guide section for error code explanations.

TROUBLESHOOTING GUIDE

In the event of an error code, unplug the 3M Attest Auto-reader 390 and plug it back in to initiate the self-diagnostic test. If error code re-occurs, call for repair or replacement.

The LCD display panel will indicate detected malfunctions by displaying a caution or error code under affected BI wells. The audible alarm, if active, sounds and the LCD display backlight flashes to alert the user to all caution (except C1) and error codes. The alarm and backlight can be silenced and extinguished, respectively, by pressing the Buzzer Alarm OFF button . Note that disabling the audible alarm does not disable the flashing LCD display backlight. Refer to the following:

Code/ Issue	Problem	Solution
C1	Unit Warm-up Caution Incubation block not at temperature.	Allow a 30 minute warm-up after 3M Attest Auto-reader 390 is plugged in. (Note: 3M Attest Rapid Readout BI readings are disabled until incubator block reaches 60°C +/- 2°C).
C2	BI Removed from Well Caution 3M Attest Rapid Readout BI has been removed before incubation is complete.	Replace 3M Attest Rapid Readout BI in proper well within 10 seconds to prevent loss of collected data.
C4	Ambient Light Caution	Move 3M Attest Auto-reader 390 to reduced ambient light conditions.
E1	Temperature Control Error This error will occur if system is no	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or

Code/ Issue	Problem	Solution
	longer able to control the temperature of the incubator block (60°C +/- 2°C).	replacement.
E3	Memory Error 3M Attest Auto-reader 390 unable to save/read data from memory.	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement.
E7	Communications Error This error will occur if internal communications have failed.	Call for repair or replacement.
E8	LED Error This error occurs for a specific well if the system detects that a UV LED or Optics Detection circuit is not working.	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement. Do not use affected incubation well.
E11	BI Inserted Incorrectly Error Incubation time remaining will not display and fluorescent readings are not taken.	Ensure BI is designed to be incubated and read in the 3M Attest Auto-reader 390. Ensure BI is fully inserted into incubation well.
E12	Power Failure Error Power to Auto-reader was off for too long to provide valid result for affected BI(s).	Ensure power supply connection to Auto-reader is secure.
E14	Ambient Light Error (occupied well) Fluorescent result will not display.	Rely on optional visual pH color change result for BI in affected well. Move 3M Attest Auto-reader 390 to reduced ambient light conditions to avoid future occurrences.
Other	Unwanted alarm sounds when positive BI detected.	Disable alarm; refer to 3M Attest Auto-reader 390 Configuration section for instructions.
Other	Alarm Sounds 3M Attest Rapid Readout Biological Indicator has been removed from incubation/reader well prematurely.	Replace 3M Attest Rapid Readout Biological Indicator in incubation/reader well within 10 seconds to prevent loss of data collected.
Other	LCD panel does not display properly.	Unplug the unit and plug it back in. If problem persists call for repair or replacement.
Other	Negative Control Vial BI not properly activated/ media did not contact spore strip.	Ensure BI cap is closed, and glass ampoule is crushed. Tap BI until media flows to and wets spore strip at bottom of BI.
Other	Negative Control Vial Control vial was sterilized.	Check process indicator on BI label. Use a non-sterilized vial as a control.
Other	Unexplained Positive BI Adhesive residue on vial.	Do not place tape or labels on vial.
Other	Unexplained Positive BI Test vial is absorbing fluorescent residue from a chemical indicator or tape.	Place the vial so it does not come in direct contact with chemical indicators or tape.
Other	3M Attest Auto-reader 390 will not enter Configuration mode. Pushbutton Configuration may be disabled on 3M Attest Auto-reader 390 Web App.	On 3M Attest Auto-reader 390 Web App Configuration Page, go to Pushbutton Configuration and select the Enable radio button. Save changes to begin using new setting.
Other	3M Attest Auto-reader 390 LCD display flashes and asterisks (***) and/or E11 and/or C4 codes appear on display below unoccupied incubation wells.	Move 3M Attest Auto-reader 390 away from sunlight or sources of incandescent light.

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Connecting to the 3M Attest Auto-reader 390 Web App

To use the 3M Attest Auto-reader 390 Web App, the 3M Attest 390 Auto-reader must be connected to the facility network using the supplied Ethernet cable. Consult with your facility IT personnel for assistance with network connectivity. Users may connect to the 3M Attest Auto-reader 390 Web App using a URL based upon either the 3M Attest Auto-reader 390 serial number or the IP Address. To view the 3M Attest Auto-reader 390 Web App from a wireless device such as a mobile hand-held device or a computer on a sub-net different from that to which the 3M Attest Auto-reader 390 is connected, use the IP Address URL. **Note:** To access the 3M Attest Auto-reader 390 Web App within a secure network, the PC or mobile device must have access to the secure network; options include connecting via a network Ethernet cable, signing in to the facility's wireless network, or connecting via a virtual private network (VPN).

A. Serial number URL

1. Connect the supplied Ethernet cable from the 3M Attest Auto-reader 390 to an Ethernet jack.
2. Depress the  button once to display the 3M Attest Auto-reader 390 6-digit serial number on the left side of the top row of text on the LCD display.
3. Enter <http://attestxxxxxx> where xxxxxx represents the unit's serial number (e.g., <http://attest001127>) from Step 2 in PC web browser address bar and press Go.
4. 3M Attest Auto-reader 390 Web App Status page will display on PC. For future ease of use, add this page to your web browser "Favorites."

B. IP Address URL

The 3M Attest Auto-reader 390 Web App may also be accessed using either a dynamic or static IP address.

1. Dynamic IP Address

A dynamic IP address is assigned when the Ethernet cable is connected from the 3M Attest Auto-reader 390 to a live Ethernet jack.

Press and release the  button twice to display the 3M Attest Auto-reader 390's dynamic IP address. Enter <http://IPaddress> (e.g. <http://123.45.678.90>) in the web browser address bar of the PC or mobile hand-held device and press Enter. The 3M Attest Auto-reader 390 Web App Status page will display.

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
You may wish to label the Ethernet jack with the assigned address. Enter <http://StaticIPaddress> (e.g. <http://111.22.333.44>) in the web browser address bar on your PC or mobile hand-held device and press Enter. The 3M Attest Auto-reader 390 Web App Status page will display.

About the 3M Attest Auto-reader 390 Web App Pages

The 3M Attest Auto-reader 390 Web App includes four pages:

a) Status page

1. User-entered nickname (optional)
2. Well Configuration (1-5 and 6-10)
3. Well Number
4. BI Result/Remaining Incubation Time
5. Date and Time
6. 3M Attest Auto-reader 390 Serial Number

The Status page provides a real-time representation of the 3M Attest Auto-reader 390 configuration and LCD panel display. The current date and time and 3M Attest Auto-reader 390 serial number are also indicated on this page. From this page, the user may navigate to:

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Note 1: The top four rows (highlighted in green) are populated automatically by the 3M Attest Auto-reader 390 and cannot be changed by the user. The remaining rows are completed by the user to satisfy the documentation requirements of their facility.

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The 3M Attest Auto-reader 390 stores results for the last 100 incubated BIs on embedded software. The 3M Attest Auto-reader 390 Web App History page provides a summary of these last 100 results.

d) Configuration page

This page allows the configuration of the 3M Attest Auto-reader 390's network and language settings from the 3M Attest Auto-reader 390 Web App.

Documenting BI Results using the 3M Attest Auto-reader 390 Web App

3M										
3M™ Attest™ Auto-reader 390										
<u>Status</u> <u>History</u>		ATTEST 1291					ATTEST 1292			
1	2	3	4	5	6	7	8	9	10	
+	59	-	-		+	-	-	174		
 2012-10-18 15:27  200004										

1. Open the Results Detail page

After inserting a new 3M Attest Rapid Readout BI in the 3M Attest Auto-reader 390, it is a best practice to open the 3M Attest Auto-reader 390 Web App Results Detail page for that BI by clicking on the appropriate well number from the Status page.

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Embedded software in the 3M Attest Auto-reader 390 saves the data from the auto-populated (i.e., green highlighted) fields. The application relies on the temporary storage in the user's web browser to save the data in the user-populated fields. In both cases, data is only stored for the last 100 records, i.e., the 101st incubated BI will erase the information saved for the 1st incubated BI. To retain the Results Detail information for future reference it is best practice to print the page or save to an electronic file as soon as the BI result is available.

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Using the Configuration Page

Certain set-up and configuration options are available on the password-protected 3M Attest Auto-reader 390 Web App Configuration page. Ideally, the 3M Attest Auto-reader 390 should be empty (i.e., have no incubating BIs) when it is reconfigured.

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1) <http://attestxxxxxx/protect/config.htm> (where xxxxxx is the 6-digit unit serial number), or
 2) <http://xxx.xx.xxx.xx/protect/config.htm> (where xxx.xx.xxx.xx is the IP address)
 in the web browser address bar and pressing Go or Enter. A security screen will appear. Fill in the fields as follows:

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Auto-reader Type	Select the configuration desired for wells 1-5/6-10 by clicking on the appropriate radio button.
Language	Selected setting indicates language in which text on 3M Attest Auto-reader 390 Web App pages will appear on the PC/mobile device screen.
Pushbutton Configuration	Selecting Enable ☐ allows the user to configure the 3M Attest Auto-reader 390 from the buttons on the unit itself. Selecting Disable ☐ prevents the user from being able to configure the 3M Attest Auto-reader 390 from the buttons on the unit itself.
Buzzer Configuration	Selecting Enable ☐ enables the audible alarm. The speaker icon will appear on the 3M Attest Auto-reader 390 LCD panel. Selecting Disable ☐ disables the audible alarm. The speaker icon will disappear from the 3M Attest Auto-reader 390 LCD panel.
Attest Auto-reader 390 Nickname	Naming the 3M Attest Auto-reader 390 may be useful in facilities with multiple units. The assigned nickname will appear on the 3M Attest Auto-reader 390 Web App Status page.
Enable DHCP	Default: Setting is checked indicating the Web App will use the Dynamic IP address automatically assigned when the 3M Attest Auto-reader 390 is connected to the network. Remove the check if you wish to use a Static IP address assigned by your facility IT.
IP Address	Enter Static IP address assigned by IT.
Gateway	Enter Gateway.
Subnet Mask	Enter Subnet Mask.
Primary DNS	Enter Primary DNS.
Use Computer Time	Check this box if you wish to set the clock on-board the 3M Attest Auto-reader 390 using the computer's clock.
Attest Auto-reader 390 Time	Alternatively, leave the Use Computer Time box unchecked and complete the date and time boxes to set the 3M Attest Auto-reader 390 clock.

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TROUBLESHOOTING THE 3M ATTEST AUTO-READER 390 WEB APP

Issue	Problem	Solution
Unable to access 3M Attest Auto-reader 390 Web App on networked computer	3M Attest Auto-reader 390 and/or computer not connected to network.	Ensure 3M Attest Auto-reader 390 Ethernet cable is connected to a live jack and that connection to 3M Attest Auto-reader 390 is secure. Ensure computer is connected to network. If failure is with connection attempted using 3M Attest Auto-reader 390 serial number URL, try connecting using IP address URL. Contact facility IT or 3M if problem persists.
Unable to access 3M Attest Auto-reader 390 Web App on mobile device	3M Attest Auto-reader 390 not connected to network. For security reasons, facility firewall may prevent mobile devices not connected to the	Ensure 3M Attest Auto-reader 390 Ethernet cable is connected to a live jack and that connection to 3M Attest Auto-reader 390 is secure. Contact facility IT for their recommendation on connecting to the secure network. Solutions may include enabling WiFi on the mobile device and signing into the facility network or connecting to

	facility network from connecting to secure IP addresses.	the facility virtual private network (VPN) on the mobile device.
Data entered by user on Results Detail page disappears	User entered data is written to temporary web browser storage. Browser may have cleared data.	Set web browser to accept (allow or enable) cookies. As a best practice, print or save Results Detail record as soon as BI result is final.
BI result displayed on 3M Attest Auto-reader 390 LCD panel does not appear on Results Detail page	The Results Detail page has not refreshed yet.	Use your browser to refresh the Results Detail page, or wait up to 60 seconds for the page to automatically refresh.
Results Detail page will not open from the Status page	The Results Detail page only opens from the Status page for occupied 3M Attest Auto-reader 390 wells.	Do not attempt to navigate to the Results Detail page from unoccupied wells on the Status page.

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 - The WEEE logo on the packaging, the appliance and in the instructions indicates that disposal must be carried out in accordance with the EU Electronic and Electrical Equipment Directive.

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2. has been signed by (a été signé par) **Ryan J Bartz**

3. acting in the capacity of (agissant en qualité de) **Notary Public, State of Minnesota**

4. bears the seal / stamp of (est revêtu du sceau / timbre de)

Ryan J Bartz, Notary Public, State of Minnesota

Certified Attesté

5. at (à) **St. Paul, Minnesota**

6. Date (Date) **02/27/2020**

7. by (par) **Secretary of State, State of Minnesota**

8. File No (Sous n°) **1145055300035**

9. Seal / Stamp (Sceau / Timbre):

10. Signature (Signature):

Steve Simon

Steve Simon

Secretary of State
State of Minnesota

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.


To verify the issuance of this Apostille, see <https://apostille.sos.state.mn.us>.

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the convention. In such cases, the certificate should be presented to the consular section of the mission representing that country.

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Operator Manual for 3M™ Attest Rapid Readout EO Biological Indicators Incubator (Auto-reader) is truthful and accurate.

Regina Feferman
International Regulatory Affairs Associate
3M Health Care

2/26/20

Date

State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 26 day of February 2020.

Ryan Bartz, Notary Public
My Commission Expires: January 31, 2024

Attest Biological Indicator Incubators
Attest Rapid Readout EO Biological Indicators Incubator
(Auto-reader)

Read and understand all safety information before using this product.

Content Disclaimers

Pictorial Disclaimer

Sample printouts, graphics, displays and screens are for information and illustration purposes only, and shall not be used for clinical or maintenance evaluations. Data shown in sample printouts and screens do not reflect actual patient names or test results.

Hardware Disclaimer

The 3M Attest Auto-reader 390G hardware specifications are subject to change.

The system images, hardware components and hardware specifications included in the manual may not match the system installed. Any changes or modifications to the authorized system installation have been verified as compatible with the functionality outlined in this document.

Product name

Attest Rapid Readout EO Biological Indicators Incubator

(Auto-reader) (Hereinafter, the devices may be referred to as "Biological Indicators Incubator (etc.), "Incubator" and etc.)

Intended Use

Biological Indicator Incubators are designed to monitor the effectiveness of the process of steam and gas sterilization of medical devices in gravitational steam and gas sterilizers.

The 3M Attest Auto-reader 390G is designed to incubate and automatically read the 3M Attest Rapid Readout Biological Indicator for Ethylene Oxide, 1294, at 37°C for a final fluorescent result at 4 hours.

Contraindications

Contraindications are not available

Scope of use

Only health care professionals should use this equipment.

Possible side effects

No known possible side effects.
supplied by 3M.

Product description

The 3M Attest Auto-reader 390G (Attest Rapid Readout EO Biological Indicators Incubator) is designed to incubate (37°C) and automatically read the 3M Attest Rapid Readout Biological Indicators (RRBIs) for a final fluorescent reading. The Attest Auto-reader 390G is also designed to allow further incubation of the compatible Attest RRBI for an optional visual pH color change indication; refer to the specific RRBI labeling and instructions for use for exact details.

The function of the device is two-fold. First, it serves as an incubator for the time period required to obtain either a rapid, fluorescent readout or to observe an optional color change result. Second, the Auto-reader determines whether a potential sterilization cycle failure may have occurred by measuring the intensity of fluorescence of the BI.

It is designed to automatically read 3M Attest Rapid Readout Biological Indicators for a final negative fluorescence reading at 4 hours for 1294 biological indicators.

Serial Number

For easy identification, each 3M Attest Auto-reader 390G has a unique serial number printed on a label found on the back of the unit and displayed on the left side of the upper row of text on the LCD Display when the user depresses the  button. Please record your serial number in this manual for future reference:

SAFETY AND PRECAUTIONS

The 3M Attest Auto-reader 390G and its related devices and accessories are designed to provide safe and reliable service when used according to the instructions provided. Please read, understand, and follow all safety information contained in the instructions for use included with the 3M Attest Auto-reader 390G and 1294 Rapid Readout Biological Indicator devices prior to use. Use this equipment only for the purpose described in this Operator's Manual. **Retain these instructions for future reference.**

The unit is designed to be used only with the Rower Supply module and Ethernet cable supplied by 3M. If this product is used in a manner not specified, the protection provided by the product may be impaired.

Explanation of Single Word Consequences

 **Warning:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

 **Caution:** Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury.

Warnings and Safety Precautions

The following warnings and precautions should be followed to avoid unsafe actions that could result in personal injury or damage to the instrument.

	WARNING: To reduce the risk associated with hazardous voltage
Use indoors only.	
Do not use the equipment if it is not working properly or if it has suffered any damage.	
Use only the power supply specified for this product and certified for the country of use.	

	CAUTION: To reduce the risk of injury or instrument damage
Do not spill liquid onto or into the instrument. Do not immerse the unit in liquid.	
Always unplug the 3M Attest Auto-reader 390G and allow to cool before cleaning.	
Clean external surfaces using only the instructions provided by the manufacturer.	
Do not open the instrument housing – there are no user serviceable parts.	
The instrument must be returned to the manufacturer for repair.	
Allow the 3M Attest Rapid Readout Biological Indicator to cool for the recommended time period before activating. Activating or excessive handling of the biological indicator before cooling may cause the glass ampoule to burst.	
Wear safety glasses when activating the 3M Attest Rapid Readout Biological Indicator.	

	WARNING: To reduce the risk associated with incorrect results
Instrument is to be used by operators familiar with the device, device functionality, and Operator's Manual.	
Do not place the instrument in environment exposed to sunlight or strong incandescent light.	
Do not place instrument close to any device that emits a strong electro-magnetic field.	
Do not use on a vibrating surface.	
Match the cap color of the 3M Attest Rapid Readout Biological Indicator with the color-coded configuration sticker surrounding the 3M Attest Auto-reader 390G incubation wells.	
Do not remove or change placement of 3M Attest Rapid Readout Biological Indicator once it is placed into a well.	
Do not remove the 3M Attest Rapid Readout Biological Indicator from the incubation well until the (+)	

or (-) symbol on the LCD panel indicates the test is complete.

Caution: To reduce the risk associated with incorrect results

To avoid the possibility of the 3M Attest Rapid Readout Biological Indicator vial absorbing fluorescent residue from a chemical indicator or tape, place the 3M Attest Rapid Readout Biological Indicator vial so it does not come in direct contact with chemical indicators or tape.

HEALTH AND SAFETY INFORMATION

Device Safety Compliance

The 3M Attest Auto-reader 390G complies with the following standards as demonstrated by the CB Scheme Certificate and test report issued by Underwriters Laboratories (UL):

- IEC 61010-1 (2001) Second Edition. Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements
- IEC 61010-2-010 (2003) Second Edition. Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials

The 3M Attest Auto-reader 390G is listed as Laboratory equipment and carries the UL mark with adjacent indicators "C" and "US" based on compliance to the standards UL 61010-1 and CAN/CSA 22.2 No. 61010-1.

The 3M Attest Auto-reader 390G complies with the CE mark related to the Low Voltage Directive (LVD) 2006/95/EC as confirmed in the Declaration of Conformity.

The 3M Attest Auto-reader 390G complies with the RoHS Directive, Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

The 3M Attest Auto-reader 390G complies with the WEEE Directive, Directive 2012/19/EU of the European Parliament and of the Council of 04 July 2012 on waste electrical and electronic equipment (WEEE).

EMC Compliance

The 3M Attest Auto-reader 390G complies with the following EMC standards as confirmed in the Certificate of Compliance generated by 3M:

IEC 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

The 3M Attest Auto-reader 390G complies with the EMC requirements of the CE mark EMC Directive 2004/108/EC.

The 3M Attest Auto-reader 390G complies with the Australian and New Zealand electrical safety and electromagnetic compatibility requirements as confirmed in the Supplier's Declaration of Conformity that is linked to the Australian/New Zealand RCM (Regulatory Compliance Mark).

Note: This equipment has been tested and found to comply with the limits for a Class A digital device pursuant to Part 15, Subpart B, of the FCC Rules. These limits are designed to provide a reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their own expense.

This Class A digital device meets all requirements of the Canadian Interference-Causing Equipment Regulations.

Class 1 LED Product

An investigation was conducted on the 3M Attest Auto-reader 390G per the requirements of IEC 60825-1 Ed. 1.2 (A2:2001). The 3M Attest Auto-reader 390G is a Class 1 LED product with Class 1 Internal Radiation.

The LED apertures are located at the bottom of each Incubator/Reader well perpendicular to the opening of each well.

NOTE: Class 1 products are safe for unprotected exposure of the eyes.

Caution - Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

Note: The use of the term "controls" in the above caution statement as required by IEC 60825-1 Ed. 1.2 (A2:2001) refers to mechanisms by which the 3M Attest Auto-reader 390G may be controlled and not to the use of biological indicator controls

EXPLANATION OF PRODUCT AND PACKAGING LABEL SYMBOLS

Attention - Refer to the instructions for use

Waste Electrical and Electronic Equipment (WEEE)

UL Listed to US and Canadian Safety Standards

Mark of Conformity to European Directives

Direct current

Compliant to all applicable ACMA regulatory arrangements (RCM)

Manufacturer

Catalog number

Authorized representative for the European Community

Serial number

Date of manufacture

The number of items in the package

Class II Equipment

BSMI (Taiwan) Marking indicates compliance with mandatory safety standards and / or electromagnetic compatibility for electrical and electronic products under the Goods Inspection Act, adopted by the Bureau of Standards, Metrology and Inspection (BSMI)

Chinese CCC certification. The so-called 3c certification is the system of mandatory product certification in China, the English name "Chinese Compulsory Certification", English abbreviation - CCC

South Korea - South Korea's KC certification requires that manufacturers of electrical / electronic equipment have a KC certification mark.

PSE marking (Japan). The marking applies to all of the specified electrical appliances and materials (category A, 115 points) in the Electrical Appliance and Material Safety (DENAN) Act

The FCC Declaration of Conformity or the FCC label, as well as the Federal Commission labeling, is a certification mark used on electronic products manufactured or sold in the United States that certifies that the electromagnetic interference level of a device is within the limits approved by the Federal Communications Commission.

This indicates that the power supply meets the Level VI energy efficiency standards required by the US Department of Energy.

Listed in UL U.S. Safety Standards

Product certification to IEC-based standards in Argentina

The GS marking is a German national mark indicating that the product has been tested and found to conform to the standards for the product.

IRAM Safety Mark License - IRAM & S Mark (System 5) in Argentina


TUV – certification body. GS – German national sign, demonstrating that the device was tested and installed and is compliant with the standards for this type of device


This symbol indicates the polarity of the plug pins:

In this case, the center pin is positive.

This indicates that the power supply is for indoor use only.

The 3M Attest Auto-reader 390G is not designed to incubate 3M Attest 1264 EO Biological Indicators, which have a sterilant entry port in the center of the top of the cap. The 3M Attest Auto-reader 390G is designed to incubate and read only 3M Attest Rapid Readout Biological Indicators for EO, catalog number 1294. Other biological indicators (BIs) are not compatible with this device and cannot be used.

Packaging material description

Biological Indicator Incubator, power supply, cord, and Operator's Manual are packaged for protection during shipping. All items are packaged inside a corrugated box.

Contents of delivery

1. Biological Indicator Incubator – 1 PC.
2. Power supply – 1 PC.
3. Power adapter – 4 PC.
4. Ethernet cable – 1 PC.
5. Operator's manual – 1 PC.

SPECIFICATIONS*

* The tolerance is $\pm 5\%$ for all values, unless otherwise specified

Parameter/characteristic	390
Dimensions, mm	350 x 180 x 130
Weight, kg	No more than 2
The number of cells for incubating indicators	10
Incubation temperature	37°C $\pm$ 2
Incubation time	2 h If necessary, additional research: incubation - 7 hours, determine the color change pH indicator
Application area	Ethylene oxide sterilization
The function of self-diagnosis of the internal temperature of the incubator	+
Function of automatic reading of biological indicators (negative fluorescence)	+
Power Supply Specifications	
Input characteristics (power supply)	
Voltage Range	100-240 Volts AC

Parameter/characteristic	390
Frequency	50/60 Hertz
Current	0,6 Amps
Output characteristics (power supply)	
Voltage Range	12 Volts DC
Current	1,5 Amps
Power supply cord length	150 cm
Battery (subject to availability)	
Type	CR 2032
Max reverse charging current	10mA
Software	
Software name	390G
Carrier	3M Company
Functionality	Incubation of BI and measurement
Version	3.0.6
Release date	June 16, 2015
Classification (in accordance with IEC 62304)	Class B
Volume	The code size of the master PIC 128k bytes and the size of the RAM 3808 bytes

Environmental Operating Conditions

Environmental Condition	Operating Condition
Altitude	3000 (max) Meters
Operating Temperature	16 – 29 °C
Operating Relative Humidity	20 – 80% (non-condensing)
Voltage Range	12 Volts DC
Frequency	DC
Current	1,2 Amps
Degrees of protection provided by enclosures (IP Code)	IPX0 (no moisture protection inside the housing)
Electrical safety class	Class I
Installation/Over Voltage	Category II
Pollution Degree	2

INSTRUCTIONS FOR USE – Initial Set-up

Choose Location

Place 3M Attest Auto-reader 390G on a firm and level surface, away from sunlight and incandescent light. To save counter space, units may be stacked two high.

WARNING:

- Do not place the instrument in environment exposed to sunlight or strong incandescent light.
- Do not place instrument close to any device that emits a strong electro-magnetic field.
- Do not use on a vibrating surface.

Select & Install Correct Power Plug

The 3M Attest Auto-reader 390G power supply comes with multiple blade configurations. Install the blade configuration appropriate for your country into the power supply and plug into appropriate AC grounded outlet.

WARNING: Use only the power supply specified for this product and certified for the country of use.


Wait 30 Minutes Before Using

Once plugged in, allow a 30 minute warm-up period before placing 3M Attest Rapid Readout Biological Indicators into the incubation wells. The C1 caution code will disappear when proper incubation temperature is reached.

Special Considerations

It is recommended that the unit be left on to eliminate warm-up periods. The unit may be turned off when it is anticipated that it will not be used for an extended period.

INSTRUCTIONS FOR USE – Wells, Buttons, & Display Interpretation

Please familiarize yourself with the face of the 3M Attest Auto-reader 390G:

Wells

1. Incubation/Reader Wells – The 3M Attest Auto-reader 390G has ten incubation wells.

Buttons

2. Information / Up Arrow Button

3. Temperature / Down Arrow Button

4. Buzzer Alarm OFF / Enter Button

Information Display

5. LCD Display – Display is backlit.

6. LCD Display – Example of Error Message for well #6

7. LCD Display – Example of Remaining Incubation Time, in minutes, for well #10

8. LCD Display – Example of Negative (-) BI result for well #9

9. LCD Display – Example of Positive (+) BI result for well #8

10. LCD Display – Alarm Enabled Indicator

Color-Coded Configuration Indicator

11. Incubation well configuration sticker is color-coded green

BI Crusher

12. BI Crusher Well

Connectors (rear panel)

- Ethernet connector (left) for connecting an Ethernet cable supplied by «3M».

- Connector (right) 2.1 mm Barrel Jack for connecting the power supply supplied by «3M».

INSTRUCTIONS FOR USE – 3M Attest Auto-reader 390G Configuration

Note: Ideally, the 3M Attest Auto-reader 390G should be empty (i.e., have no incubating BIs) when it is reconfigured.

Press the and buttons on the front of the 3M Attest Auto-reader 390G simultaneously for 3 seconds

to enter configuration mode. Once in configuration mode, use the Up  and Down  buttons to set the 3M Attest Auto-reader 390G clock, audible alarm, and incubation well configuration as desired. The button  acts as an "Enter" key and moves the flashing cursor to the next changeable setting.

Setting the 3M Attest Auto-reader 390G Clock

The 3M Attest Auto-reader 390G date and time are displayed in this format:

YYYY-MM-DD HH:MM, with the time displayed using a 24 hour clock. To change the year, if necessary, push the up or down buttons to set desired year, then press  to advance to the month, adjust as necessary and press  to advance to the day. Use the  button to move to the time and set the hours and minutes.

Setting the Audible Alarm Function

Move across the top row of text by pushing the  button. The flashing "B" in Buzzer indicates this is a changeable setting. It is recommended this be set at "Buzzer Enabled" to audibly alert the operator to a positive BI but it may be set to "Buzzer Disabled" by using either the up or down button.

Note: If Buzzer Disabled is selected, and the changes accepted, the speaker icon () will disappear from the LCD panel and the alarm will be disabled for all positive 3M Attest Rapid Readout Biological Indicators, error codes, and caution code alarm warnings.

Note: The 3M Attest Auto-reader 390G display backlight flashes whenever an audible alarm sounds. Disabling the audible alarm does NOT disable the LCD display backlight flashing.

Press  to select the desired option and advance to the Accept screen.

Accept Screen

Use the up or down buttons to toggle between Yes and No, press  to accept changes (Yes selected) or retain current settings (No selected) and exit the configuration mode. The LCD panel will return to the standard display.

INSTRUCTIONS FOR USE - Readout Procedure

Instructions for Biological Indicator (BI) Activation and Incubation

 **CAUTION:** To avoid the possibility of the 3M Attest Rapid Readout Biological Indicator vial absorbing fluorescent residue from a chemical indicator or tape, place the 3M Attest Rapid Readout Biological Indicator vial so it does not come in direct contact with chemical indicators or tape.

 **WARNING:** Match the green cap color of the 3M Attest 1294 Rapid Readout Biological Indicator with the green color-coded configuration sticker surrounding the 3M Attest Auto-reader 390G incubation wells.

 **CAUTION:** Wear safety glasses when activating 3M Attest Rapid Readout Biological Indicators.


Follow this procedure for every 3M Attest 1294 Rapid Readout Biological Indicator to be read.


While wearing safety glasses, **CLOSE** the cap of the cooled 3M Attest 1294 Rapid Readout Biological Indicator (BI) by pushing cap down.

CRUSH the glass media ampoule in the Crusher Well built into the top of the 390G Auto-reader.

Holding BI by the cap, **TAP** the bottom of the vial on a tabletop until media wets spore strip at bottom of vial. Do not tap BI on Auto-reader.

PLACE the crushed BI in **color-coded** incubation well. Observe 3M Attest Auto-reader 390G LCD panel to confirm correct placement of the BI. If the BI is correctly seated, the 3M Attest Auto-reader 390G LCD panel will display ***, **, * and then 240 remaining minutes of incubation directly below the incubating BI.
RECORD BI result.

WARNING: Do not remove the 3M Attest 1294 Rapid Readout Biological Indicator from the incubation well until the (+) or (-) symbol on the LCD panel indicates the test is complete.

INSTRUCTIONS FOR USE – Additional Information

BI Placement

Remaining incubation time, in minutes, will be displayed on the LCD display to indicate the 3M Attest Rapid Readout Biological Indicator is incubating in the reader well and should not be touched. The time remaining display indicates the results are not available and incubation is in progress.

Do not remove 3M Attest Rapid Readout Biological Indicator until the (+) or (-) symbol on the LCD panel indicates the test is complete.

Note: If 3M Attest Rapid Readout Biological Indicator is inadvertently removed before the incubation period is complete, an alarm will sound and a C2 Caution code will be displayed on the LCD panel under that well. **You must return the 3M Attest Rapid Readout Biological Indicator to the reader well within 10 seconds or the results will be lost.**

WARNING: Do not remove or change placement of 3M Attest Rapid Readout Biological Indicator once it is placed into a well. Doing so may result in loss of results and the 3M Attest Rapid Readout Biological Indicator test may be invalid.

INSTRUCTIONS FOR USE – Positive Control

Positive Control

Activate and incubate at least one non-processed 3M Attest 1294 Rapid Readout Biological Indicator to use as the positive control each day that a processed 3M Attest 1294 Rapid Readout Biological Indicator is incubated and read by the unit. Write a "C" (for "control") and the date on the label. The positive control can be placed in any incubation well. The positive control should be the same part number and have the same lot number as the processed 3M Attest Rapid Readout Biological Indicator in the unit.

The positive control biological indicator helps confirm:

- correct incubation temperatures are met (heater block reading of 37°C +/-2°C);
- viability of spores has not been altered due to improper storage temperature, humidity, or proximity to chemicals;
- capability of media to promote rapid growth; and
- proper functioning of 3M Attest Auto-reader 390G components.

INSTRUCTIONS FOR USE – Interpretation of Results

Fluorescent Results

After placing a 3M Attest Rapid Readout Biological Indicator into an incubation well, the 3M Attest Auto-reader 390G LCD panel will display ***, ", * and then the remaining minutes of incubation directly below the incubating BI to indicate that incubation/read ing is in progress.

A (+) will display on the LCD panel directly under the BI and an alarm will sound (if the sound feature is enabled) as soon as a positive 3M Attest Rapid Readout Biological Indicator result is detected.

At the end of the four hour incubation time, if a negative 3M Attest Rapid Readout Biological Indicator result is detected, a (-) will display on the LCD panel directly under the incubating BI. The final negative fluorescent result (-) for a processed 3M Attest 1294 Rapid Readout Biological Indicator after 4 hours of incubation indicates an acceptable sterilization process. The positive control should read positive (+).

After either the + or - symbol is displayed, and the result is recorded, the 3M Attest Rapid Readout Biological Indicator can be removed from the well. The 3M Attest

Rapid Readout Biological Indicator should not be removed or touched when the display indicates remaining incubation time.

The positive (unprocessed) control 3M Attest 1294 Rapid Readout Biological Indicator must provide a positive fluorescent result (+). Processed 3M Attest Rapid Readout Biological Indicator results are not valid until the positive control reads fluorescent positive (+). With a processed 3M Attest Rapid Readout Biological Indicator, a positive (+) result indicates a sterilization process failure.

Act on any positive test results immediately. Determine the cause of the positive 3M Attest Rapid Readout Biological Indicator and retest the sterilizer following current facility policies and procedures.

Optional Visual pH Color Change Result

3M Attest 1294 Rapid Readout Biological Indicators may be further incubated for a visual pH color change. Due to the high sensitivity of the fluorescent result, there is no advantage to incubating the 3M Attest 1294 Rapid Readout Biological Indicators beyond 4 hours. This optional further incubation is not intended for routine use but may be used for special studies. A humidified incubator will be required to avoid media dry out. In the case of a positive control BI, a yellow color change of the growth media will appear within 7 days. Any observation of a yellow color within the vial indicates a positive result. The pH color change result can be recorded and the BI discarded as soon as the yellow color is observed and the result recorded, in the case of a processed BI, a media color change from blue-green to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media remains blue-green, can be assessed at 7 days.

FEATURES OF THE 3M ATTEST AUTO-READER 390G

Audible Alarm

An audible alarm will sound, the LCD display backlight will flash, and a (+) symbol will be displayed on the LCD panel directly under the incubating BI when a positive result is detected. An alarm will also sound if a 3M Attest Rapid Readout Biological Indicator is removed from the 3M Attest Auto-reader 390G well before a final result has been detected. Data will be lost if the 3M Attest Rapid Readout Biological Indicator is not replaced in the incubation well within 10 seconds.

Alarm On/Off

The audible alarm for positive results can be turned off for a particular 3M Attest Rapid Readout Biological Indicator. This is accomplished by pressing the "Buzzer Alarm OFF" button after the alarm has sounded. This will turn off the alarm for that BI only. Any positive 3M Attest Rapid Readout Biological Indicator results subsequently detected in other wells will result in the alarm sounding. Removing the 3M Attest Rapid Readout Biological Indicator will reactivate the alarm for that well. Pressing the "Buzzer Alarm OFF" button will not silence the alarm for error or caution code alarm warnings. To disable the alarm, see

The positive control biological indicator helps confirm:

- correct incubation temperatures are met (heater block reading of 37°C +/-2°C);
- viability of spores has not been altered due to improper storage temperature, humidity, or proximity to chemicals;
- capability of media to promote rapid growth; and
- proper functioning of 3M Attest Auto-reader 390G components.

INSTRUCTIONS FOR USE – Interpretation of Results

Fluorescent Results

After placing a 3M Attest Rapid Readout Biological Indicator into an incubation well, the 3M Attest Auto-reader 390G LCD panel will display ***, ", * and then the remaining minutes of incubation directly below the incubating BI to indicate that incubation/reading is in progress.

A (+) will display on the LCD panel directly under the BI and an alarm will sound (if the sound feature is enabled) as soon as a positive 3M Attest Rapid Readout Biological Indicator result is detected.

At the end of the four hour incubation time, if a negative 3M Attest Rapid Readout Biological Indicator result is detected, a (-) will display on the LCD panel directly under the incubating BI. The final negative fluorescent result (-) for a processed 3M Attest 1294 Rapid Readout Biological Indicator after 4 hours of incubation indicates an acceptable sterilization process. The positive control should read positive (+).

After either the + or - symbol is displayed, and the result is recorded, the 3M Attest Rapid Readout Biological Indicator can be removed from the well. The 3M Attest

Rapid Readout Biological Indicator should not be removed or touched when the display indicates remaining incubation time.

The positive (unprocessed) control 3M Attest 1294 Rapid Readout Biological Indicator must provide a positive fluorescent result (+). Processed 3M Attest Rapid Readout Biological Indicator results are not valid until the positive control reads fluorescent positive (+). With a processed 3M Attest Rapid Readout Biological Indicator, a positive (+) result indicates a sterilization process failure.

Act on any positive test results immediately. Determine the cause of the positive 3M Attest Rapid Readout Biological Indicator and retest the sterilizer following current facility policies and procedures.

Optional Visual pH Color Change Result

3M Attest 1294 Rapid Readout Biological Indicators may be further incubated for a visual pH color change. Due to the high sensitivity of the fluorescent result, there is no advantage to incubating the 3M Attest 1294 Rapid Readout Biological Indicators beyond 4 hours. This optional further incubation is not intended for routine use but may be used for special studies. A humidified incubator will be required to avoid media dry out. In the case of a positive control BI, a yellow color change of the growth media will appear within 7 days. Any observation of a yellow color within the vial indicates a positive result. The pH color change result can be recorded and the BI discarded as soon as the yellow color is observed and the result recorded, in the case of a processed BI, a media color change from blue-green to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media remains blue-green, can be assessed at 7 days.

FEATURES OF THE 3M ATTEST AUTO-READER 390G

Audible Alarm

An audible alarm will sound, the LCD display backlight will flash, and a (+) symbol will be displayed on the LCD panel directly under the incubating BI when a positive result is detected. An alarm will also sound if a 3M Attest Rapid Readout Biological Indicator is removed from the 3M Attest Auto-reader 390G well before a final result has been detected. Data will be lost if the 3M Attest Rapid Readout Biological Indicator is not replaced in the incubation well within 10 seconds.

Alarm On/Off

The audible alarm for positive results can be turned off for a particular 3M Attest Rapid Readout Biological Indicator. This is accomplished by pressing the "Buzzer Alarm OFF" button after the alarm has sounded. This will turn off the alarm for that BI only. Any positive 3M Attest Rapid Readout Biological Indicator results subsequently detected in other wells will result in the alarm sounding. Removing the 3M Attest Rapid Readout Biological Indicator will reactivate the alarm for that well. Pressing the "Buzzer Alarm OFF" button will not silence the alarm for error or caution code alarm warnings. To disable the alarm, see

configuration section.

Requirements for installation and mounting

There are no special requirements for the installation and / or installation of a medical device. Connection and use rules are described in the operation manuals.

Stackable

To save counter space, 3M Attest Auto-reader 390Gs may be stacked two high.

Information Button

When the Information button is pressed once the 3M Attest Auto-reader 390G displays two lines of text on the LCD display.

The top row indicates (left to right):

- 3M Attest Auto-reader 390G serial number
- 3M Attest Auto-reader 390G embedded software version (V X.Y.Z)
- Current date and time.

The bottom row indicates (left to right):

- The BI type configuration for wells 1 -5.
- The BI type configuration for wells 6-10.

Press the information button twice and the 3M Attest Auto-reader 390G LCD panel will display two additional lines of text.

The top row indicates:

- Dynamic IP address (if unit is connected to an Ethernet jack).

The bottom row indicates (left to right):

- Date (YYYY-MM-DD) and time that power was supplied to the 3M Attest Auto-reader 390G

Temperature Monitoring

The 3M Attest Auto-reader 390G contains an internal temperature self-diagnostic program. Should the temperature fall outside the specified range (37°C +/-2°C) an error code will be displayed. The heater block temperature may be checked by depressing the Temperature button once.

Easy Configuration

The 3M Attest Auto-reader 390G may be configured using either the buttons located on the front of the unit or using the optional 3M Attest Auto-reader 390G Web App.

Optional 3M Attest Auto-reader 390 Web App

If desired, the 3M Attest Auto-reader 390G may be connected to a PC via an Ethernet connection for additional capabilities:

- Remotely view current status of Bis in the 3M Attest Auto-reader 390G via PC or mobile hand held devices. For full functionality, the use of an HTML5 compliant browser with JavaScript enabled is recommended.
- Remotely view and configure 3M Attest Auto-reader 390G settings.
- View the last 100 BI records.
- Print out documentation using a local or networked printer and/or save documentation as an electronic file*

Note: To access the 3M Attest Auto-reader 390G Web App within a secure network, the PC or mobile device must have access to the secure network; options include connecting via a network Ethernet cable, signing in to the facility's wireless network, or connecting via a virtual private network (VPN).

* Refer to your facility's record keeping policy and procedure.

CARE AND CLEANING

CAUTION: Always unplug the 3M Attest Auto-reader 390G and allow to cool before cleaning.

Do not spill liquid onto or into the instrument. Do not immerse the unit in liquid.

Clean external surfaces using only the instructions provided by the manufacturer.

The instrument must be returned to the manufacturer for repair.

Clean the exterior surface of the 3M Attest Auto-reader 390G by wiping it with a damp cloth moistened with a mild detergent. Do not attempt to clean interior components. If additional cleaning of the wells is required, call for repair or replacement.

Note: The 3M Attest Auto-reader 390G does not contain any user serviceable parts or require any routine maintenance. Any malfunctions will be indicated by an error code. Refer to the Troubleshooting Guide section for error code explanations.

TROUBLESHOOTING GUIDE

In the event of an error code, unplug the 3M Attest Auto-reader 390G and plug it back in to initiate the self-diagnostic test. If error code re-occurs, call for repair or replacement.

The LCD display panel will indicate detected malfunctions by displaying a caution or error code under affected BI wells. The audible alarm, if active, sounds and the LCD display backlight flashes to alert the user to all caution (except C1) and error codes. The alarm and backlight can be silenced and extinguished, respectively, by pressing the Buzzer Alarm OFF button . Note that disabling the audible alarm does not disable the flashing LCD display backlight. Refer to the following

Code/ Issue	Problem	Solution
C1	Unit Warm-up Caution Incubation block not at temperature,	Allow a 30 minute warm-up after 3M Attest Auto-reader 390G is plugged in. (Note: 3M Attest Rapid Readout BI readings are disabled until incubator block reaches 37°C +/- 2°C).
C2	BI Removed from Well Caution 3M Attest Rapid Readout BI has been removed before incubation is complete,	Replace 3M Attest Rapid Readout BI in proper well within 10 seconds to prevent loss of collected data,
C4	Ambient Light Caution	Move 3M Attest Auto-reader 390G to reduced ambient light conditions.
E1	Temperature Control Error This error will occur if system is no longer able to control the temperature of the incubator block (37°C +/- 2°C).	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement.
E3	Memory Error 3M Attest Auto-reader 390G unable to save/read data from memory.	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement.
E7	Communications Error This error will occur if internal communications have failed,	Call for repair or replacement.
E8	LED Error This error occurs for a specific well if the system detects that a UV LED or Optics Detection circuit is not working,	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement. Do not use affected incubation well.
E11	BI Inserted Incorrectly Error Incubation time remaining will not display and fluorescent readings are not taken,	Ensure BI is designed to be incubated and read in the 3M Attest Auto-reader 390G. Ensure BI is fully inserted into incubation well.
E12	Power Failure Error Power to Auto-reader was off for too long to provide valid result for affected BI(s).	Ensure power supply connection to Auto-reader is secure.
E14	Ambient Light Error (occupied well) Fluorescent result will not display.	Rely on optional visual pH color change result for BI in affected well. Move 3M Attest Auto-reader 390G to reduced ambient light conditions to avoid future occurrences.
Other	Unwanted alarm sounds when positive BI	Disable alarm; refer to 3M Attest Auto-reader

Code/ Issue	Problem	Solution
	detected	390G Configuration section for instructions.
Other	Alarm Sounds 3M Attest Rapid Readout Biological Indicator has been removed from incubation/reader well prematurely.	Replace 3M Attest Rapid Readout Biological Indicator in incubation/reader well within 10 seconds to prevent loss of data collected.
Other	LCD panel does not display properly,	Unplug the unit and plug it back in. If problem persists call for repair or replacement.
Other	Negative Control Vial BI not properly activated/ media did not contact spore strip.	Ensure BI cap is closed, and glass ampoule is crushed. Tap BI until media flows to and wets spore strip at bottom of BI.
Other	Negative Control Vial Control vial was sterilized.	Check process indicator on BI label. Use a non-sterilized vial as a control.
Other	Unexplained Positive BI Adhesive residue on vial.	Do not place tape or labels on vial.
Other	Unexplained Positive BI Test vial is absorbing fluorescent residue from a chemical indicator or tape,	Place the vial so it does not come in direct contact with chemical indicators or tape.
Other	3M Attest Auto-reader 390G will not enter Configuration mode, Pushbutton Configuration may be disabled on 3M Attest Auto-reader 390G Web App.	On 3M Attest Auto-reader 390G Web App Configuration Page, go to Pushbutton Configuration and select the Enable radio button. Save changes to begin using new setting.
Other	3M Attest Auto-reader 390G LCD display flashes and asterisks (***) and/or E11 and/or C4 codes appear on display below unoccupied incubation wells.	Move 3M Attest Auto-reader 390G away from sunlight or sources of incandescent light.

OPTIONAL 3M ATTEST AUTO-READER 390G WEB APP

Web Browser Requirements

The Web Interface is designed to be viewed in Internet Explorer 8.0 or higher. Earlier versions of Internet Explorer or other browsers may result in some web pages being displayed in a less-than-optimum manner. For best results, switching to Internet Explorer 8.0 or higher is recommended. For full functionality, the use of an HTML5 compliant browser with JavaScript enabled is recommended.

Connecting to the 3M Attest Auto-reader 390G Web App

To use the 3M Attest Auto-reader 390G Web App, the 3M Attest 390G Auto-reader must be connected to the facility network using the supplied Ethernet cable. Consult with your facility IT personnel for assistance with network connectivity. Users may connect to the 3M Attest Auto-reader 390G Web App using a URL based upon either the 3M Attest Auto-reader 390G serial number or the IP Address. To view the 3M Attest Auto-reader 390G Web App from a wireless device such as a mobile hand-held device or a computer on a sub-net different from that to which the 3M Attest Auto-reader 390G is connected, use the IP Address URL. **Note:** To access the 3M Attest Auto-reader 390G Web App within a secure network, the PC or mobile device must have access to the secure network; options include connecting via a network Ethernet cable, signing in to the facility's wireless network, or connecting via a virtual private network (VPN).

A. Serial number URL

1. Connect the supplied Ethernet cable from the 3M Attest Auto-reader 390G to an Ethernet jack.
2. Depress the O button once to display the 3M Attest Auto-reader 390G 6-digit serial number on the left side of the top row of text on the LCD display.
3. Enter <http://attestxxxxxx> where xxxxxx represents the unit's serial number (e.g., <http://attest401008>) from Step 2 in PC web browser address bar and press Go.
4. 3M Attest Auto-reader 390G Web App Status page will display on PC. For future ease of use, add this page to your web browser "Favorites."

B. IP Address URL

The 3M Attest Auto-reader 390G Web App may also be accessed using either a dynamic or static IP address.

1. Dynamic IP Address

A dynamic IP address is assigned when the Ethernet cable is connected from the 3M Attest Auto-reader 390G to a live Ethernet jack.

Press and release the  button twice to display the 3M Attest Auto-reader 390G's dynamic IP address. Enter `http://IPaddress` (e.g. `http://123.46.678.90`) in the web browser address bar of the PC or mobile hand-held device and press Enter. The 3M Attest Auto-reader 390G Web App Status page will display.

2. Static IP address

Some facility IT personnel may require devices connecting to the network to use a Static IP address. If this is the case at your facility, contact the IT department to request a Static IP address.


Record the static IP address provided by your facility IT department here: _____.

You may wish to label the Ethernet jack with the assigned address. Enter `http://StaticIPaddress` (e.g. `http://111.22.333.44`) in the web browser address bar on your PC or mobile hand-held device and press Enter. The 3M Attest Auto-reader 390G Web App Status page will display.

About the 3M Attest Auto-reader 390G Web App Pages

The 3M Attest Auto-reader 390G WebApp includes four pages:

a) Status page

1. User-entered nickname (optional)
2. Well Configuration (1-5 and 6-10)
3. Well Number
4. BI Result/Remaining Incubation Time
5. Date and Time
6. 3M Attest Auto-reader 390G Serial Number

The Status page provides a real-time representation of the 3M Attest Auto-reader 390G LCD panel display. The current date and time and 3M Attest Auto-reader 390G serial number are also indicated on this page. From this page, the user may navigate to:

- the History page by clicking on History on the left of the screen.
- the Results Detail page by clicking on the well number of an occupied well (i.e., a well that currently has an incubating BI). Clicking on an unoccupied well takes the user to the History page.

b) Results Detail page

3M™ Attest™ Auto-reader 390G Results Detail

[Status](#)
[History](#)
[Print](#)

Result	
Start Time	2012-10-29 11:08:04
Well# and Configuration	#3 1294
390G Serial Number	100122
Control?	
BI Lot#	
Load#	
Sterilizer#	
Technician	

This page provides an individual record for an incubated BI. For complete documentation, the user may print a hard copy and/or save the record as an electronic file.

Note 1: The top four rows (highlighted in green) are populated automatically by the 3M Attest Auto-reader 390G and cannot be changed by the user. The remaining rows are completed by the user to satisfy the documentation requirements of their facility.

Note 2: The Results Detail page refreshes less frequently than the Status page. Until the BI result is available, the Results Detail page displays remaining minutes of incubation in the results field.

c) History page

3M™ Attest™ Auto-reader 390G

[Status](#)
[History](#)

ID	Well# and Configuration	Start Time	Result
1	#4 1294	2012-10-30 12:31:25	---
2	#3 1294	2012-10-30 09:57:01	---
3	#3 1294	2012-10-30 07:51:51	+
4	#4 1294	2012-10-30 16:15:55	-
5	#5 1294	2012-10-30 16:01:07	-
6	#1 1294	2012-10-30 15:58:18	+
7	#4 1294	2012-10-30 11:52:27	-
8	#4 1294	2012-10-30 11:42:26	C2

The 3M Attest Auto-reader 390G stores results for the last 100 incubated Bis on embedded software. The 3M Attest Auto-reader 390G Web App History page provides a summary of these last 100 results.

d) Configuration page

This page allows the configuration of the 3M Attest Auto-reader 390G's network and language settings from the 3M Attest Auto-reader 390G Web App.

Documenting BI Results using the 3M Attest Auto-reader 390G WebApp

3M

3M™ Attest™ Auto-reader 390G

Status
History

ATTEST 1294

1	2	3	4	5	6	7	8	9	10
229	-	-	-	231					

2012-10-30 16:06

SN 100122

1. Open the Results Detail page

After inserting a new 3M Attest Rapid Readout BI in the 3M Attest Auto-reader 390G, it is a best practice to open the 3M Attest Auto-reader 390G Web App Results Detail page for that BI by clicking on the appropriate well number from the Status page.

2. Enter the specific information for that BI (e.g., BI Lot #, sterilizer/load number, Technician name, etc.) in the provided editable fields. Use the Notes section to capture information not prescribed in the provided row titles; for example, you may wish to use this field to record load items.

3M

3M™ Attest™ Auto-reader 390G Results Detail

Status
History
Print

Result	223
Start Time	2012-10-30 15:58:18
Well# and Configuration	#1 1294
390G Serial Number	100122
Control?	No
BI Lot#	2014-11AD
Load#	1
Sterilizer#	2
Technician	sdf

3. The entered information is **saved** automatically to temporary web browser storage on the user's PC every 60 seconds, or by leaving the Results Details page.

4. **Print/File.** Once the BI result is available, print or file the Results Detail page to create a record. The page may be physically printed and/or saved as an electronic file.

Important! Data is only stored for the last 100 incubated BIs.

Embedded software in the 3M Attest Auto-reader 390G saves the data from the auto-populated (i.e., green highlighted) fields. The application relies on the temporary storage in the user's web browser to save the data in the user-populated fields. In both cases, data is only stored for the last 100 records, i.e., the 101st incubated BI will erase the information saved for the 1st incubated BI. To retain the Results Detail information for future reference it is best practice to print the page or save to an electronic file as soon as the BI result is available.

Using the BI History Page

Status
History

ID	Well# and Configuration	Start Time	Result
1	#4 1294	2012-10-30 12:31:25	---
2	#3 1294	2012-10-30 09:57:01	---
3	#3 1294	2012-10-30 07:51:51	+
4	#4 1294	2012-10-30 16:15:55	-
5	#5 1294	2012-10-30 16:01:07	-
6	#1 1294	2012-10-30 15:58:18	+
7	#4 1294	2012-10-30 11:52:27	-
8	#4 1294	2012-10-30 11:42:26	C2

To view or print an individual record, click on the ID number of the desired record to navigate to the Results Detail page for that BI.

Note 1: Column 1, "ID" is the assigned ID number with ID 1 representing the newest record and ID 100 representing the oldest record. For convenience, all columns are sortable; to help find a particular record, simply click once on the column heading to sort in ascending order and click a second time to reverse the order.

Note 2: The History page refreshes less frequently than the Status page. Until the BI result is available, the History page displays *** in the Results column indicating incubation is still in-process.

Using the Configuration Page

Certain set-up and configuration options are available on the password-protected 3M Attest Auto-reader 390G Web App Configuration page. Ideally, the 3M Attest Auto-reader 390G should be empty (i.e., have no incubating Bis) when it is reconfigured. Navigate to the Configuration page by entering either.

Navigate to the Configuration page by entering either:

- 1) <http://attestxxxxxx/protect/config.htm> (where xxxxxx is the 6-digit unit serial number), or
- 2) <http://xxx.xx.xxx.xx/protect/config.htm> (where xxx.xx.xxx.xx is the IP address)

in the web browser address bar and pressing Go or Enter. A security screen will appear. Fill in the fields as follows:

User name: manager

Password: Attest

and press OK.

Click on the desired radio buttons to change settings.

Language	Selected setting indicates language in which text on 3M Attest Auto-reader 390G Web App pages will appear on the PC/mobile device screen.
Pushbutton Configuration	Selecting Enable ☐ allows the user to configure the 3M Attest Auto-reader 390G from the buttons on the unit itself. Selecting Disable ☐ prevents the user from being able to configure the 3M Attest Auto-reader 390G from the buttons on the unit itself.
Buzzer Configuration	Selecting Enable ☐ enables the audible alarm. The speaker icon will appear on the 3M Attest Auto-reader 390G LCD panel. Selecting Disable ☐ disables the audible alarm. The speaker icon will disappear from the 3M Attest Auto-reader 390G LCD panel.
Attest Auto-reader 390G Nickname	Naming the 3M Attest Auto-reader 390G may be useful in facilities with multiple units. The assigned nickname will appear on the 3M Attest Auto-reader 390G Web App Status page.
Enable DHCP	Default: Setting is checked indicating the Web App will use the Dynamic IP address automatically assigned when the 3M Attest Auto-reader 390G is connected to the network. Remove the check if you wish to use a Static IP address assigned by your facility IT.
IP Address	Enter Static IP address assigned by IT.
Gateway	Enter Gateway.

Subnet Mask	Enter Subnet Mask.
Primary DNS	Enter Primary DNS.
Use Computer Time	Check this box if you wish to set the clock on-board the 3M Attest Auto-reader 390G using the computer's clock.
Attest Auto-reader 390G Time	Alternatively, leave the Use Computer Time box unchecked and complete the and time boxes to set the 3M Attest Auto-reader 390G clock.

After changing desired settings, click on Save Config button to save your changes. A new screen confirms the settings were successfully changed appears.

TROUBLESHOOTING THE 3M ATTEST AUTO-READER 390G WEB APP

Issue	Problem	Solution
Unable to access 3M Attest Auto-reader 390G Web App on networked computer	3M Attest Auto-reader 390G and/or computer not connected to network.	Ensure 3M Attest Auto-reader 390G Ethernet cable is connected to a live jack and that connection to 3M Attest Auto-reader 390G is secure. Ensure computer is connected to network. If failure is with connection attempted using 3M Attest Auto-reader 390G serial number URL, try connecting using IP address URL. Contact facility IT or 3M if problem persists.
Unable to access 3M Attest Auto-reader 390G Web App on mobile device	3M Attest Auto-reader 390G not connected to network. For security reasons, facility firewall may prevent mobile devices not connected to the facility network from connecting to secure IP addresses.	Ensure 3M Attest Auto-reader 390G Ethernet cable is connected to a live jack and that connection to 3M Attest Auto-reader 390G is secure. Contact facility IT for their recommendation on connecting to the secure network. Solutions may include enabling Wi-Fi on the mobile device and signing into the facility network or connecting to the facility virtual private network (VPN) on the mobile device.
Data entered by user on Results Detail page disappears	User entered data is written to temporary web browser storage. Browser may have cleared data.	Set web browser to accept (allow or enable) cookies. As a best practice, print or save Results Detail record as soon as BI result is final.
BI result displayed on 3M Attest Auto-reader 390G LCD panel does not appear on Results Detail page	The Results Detail page has not refreshed yet.	Use your browser to refresh the Results Detail page, or wait up to 60 seconds for the page to automatically refresh.
Results Detail page will not open from the Status page	The Results Detail page only opens from the Status page for occupied 3M Attest Auto-reader 390G wells.	Do not attempt to navigate to the Results Detail page from unoccupied wells on the Status page.

Storage and Transportation: Temperature from -10°C to +50°C, Relative Humidity Limitation from 20% to 80% (non-condensing), a pressure are not regulated.

During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

Maintenance and repair

Repairs must be done by a qualified service person. Attest incubators that require service or repair should be returned, contact the local subsidiary. The 3M Attest Incubator is serviceable electronics that does not have a designated shelf life due to the ability to repair and replace parts.

The duration of the warranty period is 12 months.

Disposal

Device contains a lithium battery. Dispose of content/container in accordance with the local/regional/national/international/regulations.

The class of medical waste is Class A (when disposing of a product that has not been used for its intended purpose).

The class of medical waste is Class B (when disposing of the product used for its intended purpose).

■ - The WEEE logo on the packaging, the appliance and in the instructions indicates that disposal must be carried out in accordance with the EU Electronic and Electrical Equipment Directive.

The product, as well as used batteries and other accessories, must not be disposed of with unsorted household waste.

Dispose of the product separately at special collection points for electronic and electrical equipment, in accordance with local legal requirements for the disposal of electronic devices.

Warranty

The 3M Attest Steam Incubator has a one year limited product warranty. The US warranty, remedy and limitations are described in the warranty card provided with the device. For countries outside of the US, any warranty is established by the subsidiary with responsibility for service of the incubator.

Legal Manufacturer

3M Company, 3M Health Care

3M Center, 2510 Conway Avenue, Bldg. 275-5W-06, St. Paul, Minnesota, 55144, USA

Authorized Representative of Manufacturer in Russia

AO «3M Russia»

108811, Moscow, p. Moskovsky, Kievskoe highway, 22nd km, house. 6, b. 1

Telephone numbers: +7 (495) 784 7474 / +7 (495) 784 7475 (fax)

e-mail: 3mra@mmm.com

Issue Date: 2019-06

Applicable national standards: A list of applicable national standards is available on request.

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State of Minnesota - Secretary of State

This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

Apostille

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**1. Country (Pays): United States of America**

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2. has been signed by (a été signé par) Ryan J Bartz**3. acting in the capacity of (agissant en qualité de) Notary Public, State of Minnesota**

**4. bears the seal / stamp of (est revêtu du sceau / timbre de)
Ryan J Bartz, Notary Public, State of Minnesota**

**Certified
Attesté**

5. at (à) St. Paul, Minnesota**6. Date (Date) 02/27/2020****7. by (par) Secretary of State, State of Minnesota****8. File No (Sous n°) 1145055300051****9. Seal / Stamp (Sceau / Timbre):****10. Signature (Signature):**

Steve Simon

Steve Simon
Secretary of State
State of Minnesota

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
3M

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Operator Manual for 3M™ Attest Rapid Readout vaporized hydrogen peroxide Biological Indicators Incubator (Auto-reader) is truthful and accurate.

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

3M Attest Biological Indicator Incubators

Attest Rapid Readout vaporized hydrogen peroxide Biological Indicators Incubator (Auto-reader)

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Hardware Disclaimer

The 3M Attest Auto-reader 490 hardware specifications are subject to change.

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Product name

Attest Rapid Readout vaporized hydrogen peroxide Biological Indicators Incubator (Auto-reader)

(Hereinafter, the devices may be referred to as "Biological Indicators Incubator (etc.)", "Incubator" and etc.)

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The 3M Attest Super Rapid Readout Biological Indicator System provides rapid, reliable monitoring of steam sterilization processes performed in health care facilities.

The system is comprised of biological indicators, challenge packs, and incubator/readers (3M Attest Auto-reader 490 and 3M Attest Auto-reader 490H having software version 4.0.0 or greater) with embedded software.

The 3M Attest Auto-reader 490 having software version 4.0.0 or greater is designed to incubate and automatically read 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V at 60°C for a final negative fluorescence reading at 24 minutes. The 3M Attest Auto-reader 490 is also designed to allow for further incubation of 3M Attest Super Rapid Readout Biological Indicators, catalog numbers 1491 and 1492V, for a final negative, visual pH color change result at 24 hours and 48 hours, respectively. This further incubation is optional and not intended for routine use.

3M Attest Auto-readers 490 having software version 4.0.0 or greater may also be used as part of the 3M Attest Rapid Readout Biological Indicator System for Vaporized Hydrogen Peroxide Sterilization to provide rapid and reliable routine monitoring of vaporized hydrogen peroxide sterilization processes in health care facilities. This system is comprised of incubator/readers (3M Attest Auto-reader 490 having software version 4.0.0 or greater and 3M Attest Auto-reader 490H) with embedded software, and a biological indicator, the 3M Attest Rapid Readout Biological Indicator 1295. The 3M Attest Auto-reader 490 having software version 4.0.0 or greater is designed to incubate and automatically read 3M Attest Rapid Readout Biological Indicators 1295 at 60°C for a final negative fluorescence reading at 24 minutes.

A positive fluorescence reading or a visual pH color change indicate a sterilization process failure. However, due to the high sensitivity of the fluorescent result, there is no advantage to continuing to incubate processed 3M Attest Biological Indicators once the fluorescent result has been documented.

Other biological indicators (BIs) are not compatible with this device and cannot be used.

The 3M Attest Auto reader 490 is designed to incubate and automatically read 3M Attest Rapid Readout Biological Indicators 1295 and 3M Attest Super Rapid Readout Biological Indicators, catalog numbers 1491 and 1492V, at 60°C for a final fluorescent result at 24 minutes.

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No known possible side effects.

Serial Number

For easy identification, each 3M Attest Auto-reader 490 has a unique serial number printed on a label found on the back of the unit and displayed on the left side of the upper row of text on the LCD Display when the user depresses the  button. Please record your serial number in this manual for future reference:

SAFETY AND PRECAUTIONS

The 3M Attest Auto-reader 490 and its related devices and accessories are designed to provide safe and reliable service when used according to the instructions provided. Please read, understand, and follow all safety information contained in the instructions for use included with the 3M Attest Auto-reader 490, the 3M Attest Rapid Readout Biological Indicator 1295, and 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V devices prior to use. Use this equipment only for the purpose described in the Operators Manual. **Retain these instructions for future reference.**

The unit is designed to be used only with the Power Supply module and Ethernet cable supplied by 3M. If this product is used in a manner not specified, the protection provided by the product may be impaired.

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Warning: Indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

Caution: Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury.

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WARNING: To reduce the risk associated with hazardous voltage

Use indoors only.

Do not use the equipment if it is not working properly or if it has suffered any damage.

Use only the power supply specified for this product and certified for the country of use.

CAUTION: To reduce the risk of injury or instrument damage

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Always unplug the 3M Attest Auto-reader 490 and allow to cool before cleaning.

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Do not open the instrument housing – there are no user serviceable parts. The instrument must be returned to the manufacturer for repair.

Wear safety glasses and gloves when activating the 3M Attest Rapid Readout Biological Indicator 1295

Allow the 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V to cool for the recommended time period before activating. Activating or excessive handling of the biological indicator before cooling may cause the glass ampoule to burst.

Wear safety glasses when activating the 3M Attest Super Rapid Readout Biological Indicators 1491 or 1492V.

WARNING: To reduce the risk associated with incorrect results

Instrument is to be used by operators familiar with the device, device functionality, and Operator's Manual.
Do not place the instrument in environment exposed to sunlight or strong incandescent light.
Do not place instrument close to any device that emits a strong electro-magnetic field.
Do not use on a vibrating surface.
Activate and incubate 3M Attest Rapid Readout Biological Indicators 1295 within 1 hour of the completion of the sterilization cycle.
Do not incubate a 3M Attest Biological Indicator if, after processing and before BI activation, it is observed to have a broken media ampoule. Retest the sterilizer with a new biological indicator.
Allow 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V to cool for the recommended time period before activating.
Do not remove or change placement of 3M Attest Biological Indicator once it is placed into a well.
Do not remove the 3M Attest Biological Indicator from the incubation well until the (+) or (-) symbol on the LCD panel indicates the test is complete.

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To avoid the possibility of the 3M Attest Biological Indicator vial absorbing fluorescent residue from a chemical indicator or tape, place the 3M Attest Biological Indicator vial so it does not come in direct contact with chemical indicators or tape.

HEALTH AND SAFETY INFORMATION

Device Safety Compliance

The 3M Attest Auto-reader 490 complies with the following standards as demonstrated by the CB Scheme Certificate and test report issued by Underwriters Laboratories (UL):

- IEC 61010-1 (2010) 3rd Edition. Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements
- IEC 61010-2-010:2014 3rd Edition. Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials

The 3M Attest Auto-reader 490 is listed as Laboratory equipment and carries the UL mark with adjacent indicators "C" and "US" based on compliance to the standards UL 61010-1, CAN/CSA 22.2 No. 61010-1 and CAN/CSA 22.2 No. 61010-2-010.

The 3M Attest Auto-reader 490 complies with the CE mark related to the Low Voltage Directive (LVD) 2006/95/EC as confirmed in the Declaration of Conformity.

The 3M Attest Auto-reader 490 complies with the RoHS Directive, Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

The 3M Attest Auto-reader 490 complies with the WEEE Directive, Directive 2012/19/EU of the European Parliament and of the Council of 04 July 2012 on waste electrical and electronic equipment (WEEE).

EMC Compliance

The 3M Attest Auto-reader 490 complies with the following EMC standards as confirmed in the Certificate of Compliance generated by 3M:

- IEC 61326-1:2012/EN 61326:2013 Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements
- EMC requirements of the CE mark EMC Directive 2004/108/EC

The 3M Attest Auto-reader 490 complies with the Australian and New Zealand electrical safety and electromagnetic compatibility requirements as confirmed in the Supplier's Declaration of Conformity that is linked to the Australian/New Zealand RCM (Regulatory Compliance Mark).

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15, Subpart B, of the FCC Rules. These limits are designed to provide a reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required

to correct the interference at their own expense.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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The FCC Declaration of Conformity or the FCC label, as well as the Federal Commission labeling, is a certification mark used on electronic products manufactured or sold in the United States that certifies that the electromagnetic interference level of a device is within the limits approved by the Federal Communications Commission. This symbol is required for China RoHS II compliance. It indicates an Environment Friendly Use Period (EFUP) of 10 years. This is the period of time before any of the RoHS substances are likely to leak out, causing possible harm to health and the environment. Each product that contains RoHS substances above the maximum permitted concentration values is labelled with a circle composed of two arrows containing a number that gives the EFUP in years.

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Contents of delivery

1. Biological Indicator Incubator – 1 PC.
2. Power supply – 1 PC.
3. Power cord – 1 PC.
4. Ethernet cable – 1 PC.
5. Biological Indicator Activator – 1 PC.
6. Operator's manual – 1 PC.

SPECIFICATIONS*

* The tolerance is $\pm 5\%$ for all values, unless otherwise specified

Parameter/characteristic	490
Dimensions, mm	228,6 x 82,55 x 222,25
Weight, kg	No more than 2
The number of cells for incubating indicators	10
Incubation temperature	60°C $\pm 2^\circ\text{C}$
Incubation time	1 h (for biological indicator 60 °C, 1491, and 1492V)
Application area	Steam (incl. vaporized hydrogen peroxide)
The function of self-diagnosis of the internal temperature of the incubator	+
Function of automatic reading of biological indicators (negative fluorescence)	+
Power cord length	230 cm
Power Supply Specifications	
Input characteristics (power supply)	
Voltage Range	100-240 Volts AC
Frequency	50/60 Hertz
Current	0,48 Amps
Output characteristics (power supply)	
Voltage Range	12 Volts DC
Current	1,6 Amps
Power supply cord length	150 cm
Voltage Range	100-240 Volts AC
Frequency	50/60 Hertz
Current	0,6 Amps
Battery (subject to availability)	

Parameter/characteristic	490
Type	CR 2032
Max reverse charging current	10mA
Software	
Software name	490
Carrier	3M Company
Functionality	Incubation of BI and measurement
Version	4.0.0
Release date	June 16, 2015
Classification (in accordance with IEC 62304)	Class B
Volume	The code size of the master PIC 128k bytes and the size of the RAM 3808 bytes

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Environmental Condition	Operating Condition	Units
Altitude	3000 (max)	Meters
Operating Temperature	16 – 40	°C
Operating Relative Humidity	20 – 80 (non-condensing)	%
Voltage Range	12	Volts DC
Frequency	DC	
Current	1,2	Amps
Degrees of protection provided by enclosures (IP Code)	IPX0 (no moisture protection inside the housing)	
Electrical safety class	Class I	
Installation/Over Voltage	Category II	
Pollution Degree	2	

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Choose Location

Place 3M Attest Auto-reader 490 on a firm level surface away from sunlight and incandescent light. To save counter space, units may be stacked two high.

Plug In Unit

Connect the power cord to the power supply. Connect the power supply to the 3M Attest Auto-reader 490 and plug the power cord into an AC grounded outlet.

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Wait 30 Minutes Before Using

Once plugged in, allow a 30 minute warm-up period before placing biological indicators into the incubation wells. The C1 caution code will disappear when proper incubation temperature is reached.

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
INSTRUCTIONS FOR USE – Wells, Buttons, & Display Interpretation

Please familiarize yourself with the face of the 3M Attest Auto-reader 490:

Wells

1. Incubation/Reader Wells – The 3M Attest Auto-reader 490 has ten incubation wells.

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10. LCD Display – Alarm Enabled Indicator

Color-Coded Configuration Indicator

11. Incubation well configuration sticker is black.

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INSTRUCTIONS FOR USE – 3M Attest Auto-reader 490 Configuration

Note: Ideally, the 3M Attest Auto-reader 490 should be empty (i.e., have no incubating BIs) when it is reconfigured.

Press the and buttons on the front of the 3M Attest Auto-reader 490 simultaneously for 3 seconds to enter configuration mode. Once in configuration mode, use the Up and Down buttons to set the 3M Attest Auto-reader 490 clock, audible alarm, and Web App as desired.

The button acts as an “Enter” key and moves the flashing cursor to the next changeable setting.

Setting the 3M Attest Auto-reader 490 Clock

The 3M Attest Auto-reader 490 date and time are displayed in this format:

YYYY-MM-DD HH:MM, with the time displayed using a 24 hour clock. To change the year, if necessary,

push the up or down buttons to set desired year, then press  to advance to the month, adjust as necessary and press  to advance to the day. Use the  button to move to the time and set the hours and minutes.

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Note: The 3M Attest Auto-reader 490 display backlight flashes whenever an audible alarm sounds. Disabling the audible alarm does NOT disable the LCD display backlight flashing.

Press  to select the desired option and advance to the next changeable setting.

Setting the Optional 3M Attest Auto-reader 490 Web App Function

It is recommended that the default **Web Pages Enabled** setting is used. If desired, however, the 3M Attest Auto-reader 490 Web App can be disabled by pressing either the up or down button. Press  to select the desired option and advance to the Accept screen.

Note: Disabling the 3M Attest Auto-reader 490 Web App will not affect connectivity to 3rd party software.

Accept Screen

Use the up or down buttons to toggle between Yes and No, press to  accept changes (Yes selected) or retain current settings (No selected) and exit the configuration mode. The LCD panel will return to the standard display.

Note: Incubation wells on 3M Attest Auto-readers 490 having software version 4.0.0 or greater are not configurable. Contact 3M for instructions on configuring incubation wells of Auto-readers having software versions less than 4.0.0.

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 **CAUTION:** Allow the 3M Attest Super Rapid Readout Biological Indicator to cool for the recommended time period before activating. Activating or excessive handling of the biological indicator before cooling may cause the glass ampoule to burst.

 **CAUTION:** To avoid the possibility of the 3M Attest Biological Indicator vial absorbing fluorescent residue from a chemical indicator or tape, place the 3M Attest Biological Indicator vial so it does not come in direct contact with chemical indicators or tape.

 **CAUTION:** Wear safety glasses when activating the 3M Attest Super Rapid Readout Biological Indicator.

 **CAUTION:** Wear safety glasses and gloves when activating the 3M Attest Rapid Readout Biological Indicator 1295.

 **WARNING:** To reduce the risk associated with incorrect results.

Do not incubate a 3M Attest Biological Indicator if, after processing and before BI activation, it is observed to have a broken media ampoule. Retest the sterilizer with a new biological indicator.

Activate and incubate 3M Attest Rapid Readout Biological Indicators 1295 within 1 hour of the completion of the sterilization cycle.


Allow 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V to cool for the recommended time period before activating.

Incubate the 3M Attest Rapid Readout Biological Indicator 1295 only in 3M Attest Auto-readers 490 having software versions 4.0.0 or greater.

Follow this procedure for every 3M Attest Biological Indicator to be read:


While wearing safety glasses (and gloves for the 1295 BI), **PLACE** BI in Attest Biological Indicator Activator.

CLOSE and **SQUEEZE** Attest Biological Indicator Activator to close BI cap and crush media ampoule.

Immediately remove BI, visually **VERIFY** the glass ampoule (containing growth media) has been crushed, and **FLICK** BI to ensure purple growth media flows to the growth chamber at the bottom of the vial.

PLACE the activated BI into any 3M Attest Auto-reader 490 well and observe 3M Attest Auto-reader 490 LCD panel to confirm correct placement of the BI. If the BI is correctly seated, the 3M Attest Auto-reader 490 LCD panel will display ***, **, * and then remaining minutes of incubation directly below the incubating BI. **RECORD** BI result.

WARNING: Do not remove the biological indicator from the incubation well until the (+) or (-) symbol on the LCD panel indicates the test is complete.

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BI Placement

All incubation wells in 3M Attest Auto-readers 490 having software version 4.0.0 or greater are programmed to provide a final fluorescent result within 24 minutes of incubation and may be used to incubate both 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V and the 3M Attest Rapid Readout Biological Indicator 1295.

Remaining incubation time, in minutes, will be displayed on the LCD display to indicate the 3M Attest Biological Indicator is incubating in the reader well and should not be touched. The time remaining display indicates the results are not available and incubation is in progress.

Do not remove 3M Attest Biological Indicators until the (+) or (-) symbol on the LCD panel indicates the test is complete. Once the result is determined, it is continuously displayed on the LCD panel until 10 seconds after the BI has been removed from the incubation well.

Note: After the first minute of incubation, if 3M Attest Biological Indicators are inadvertently removed before the incubation period is complete, an alarm will sound and a C2 Caution code will also be displayed on the LCD panel under that well. You **must return the 3M Attest Biological Indicator to the reader well within 10 seconds or the results will be lost.**

WARNING: Do not remove or change placement of 3M Attest Biological Indicator once it is placed into a well. Doing so may result in loss of results and the biological indicator test may be invalid.

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Positive Control

Activate and incubate at least one non-processed 3M Attest Biological Indicator to use as the positive control each day that a processed 3M Attest Biological Indicator is incubated and read by the unit. Write a "C" (for "control") and the date on the label. The positive control can be placed in any incubation well. A positive control should be run for each catalog number, and have the same lot number as the processed 3M Attest Biological Indicator.

It is good practice to use a positive control each day a processed 3M Attest Biological Indicator is incubated. This helps confirm:

- correct incubation temperatures are met (heater block reading of 60°C +/-2°C).
- viability of spores has not been altered due to improper storage temperature, humidity, or proximity to chemicals
- capability of media to promote rapid growth, and
- proper functioning of 3M Attest Auto-reader 490 components.

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Fluorescent Results

After placing a 3M Attest Biological Indicator into an incubation well, the 3M Attest Auto-reader 490 LCD panel will display ***, **, * and then the remaining minutes of incubation directly below the incubating BI to indicate that incubation/reading is in progress. The 3M Attest Biological Indicator should not be removed or touched when the display indicates remaining incubation time.

A (+) will display on the LCD panel directly under the BI and an alarm will sound (if the sound feature is enabled), as soon as a positive 3M Attest Biological Indicator result is detected.

At the end of the specified incubation time, if a negative 3M Attest Biological Indicator result is detected, a (-) will display on the LCD panel directly under the incubating BI. The final negative fluorescent result (-) for a processed 3M Attest Biological Indicator after 24 minutes of incubation indicates an acceptable sterilization process.

After either the + or - symbol is displayed, and the result is recorded, the 3M Attest Biological Indicator can be removed from the well. The result (+ or -) is displayed on the LCD panel until the BI is removed from the incubation well. The result clears ten seconds after the BI has been removed.

The positive (unprocessed) control 3M Attest Biological Indicator must provide a positive fluorescent result (+). Processed 3M Attest Biological Indicator results are not valid until the positive control reads fluorescent positive (+). With a processed 3M Attest Biological Indicator, a positive (+) result indicates a sterilization process failure.

Act on any positive test results immediately. Determine the cause of the positive 3M Attest Biological Indicator and retest the sterilizer following current facility policies and procedures.

Optional Visual pH Color Change Result - 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V

3M Attest Super Rapid Readout Biological Indicators may be further incubated for a visual pH color change. This optional further incubation is not intended for routine use and may be used for special studies. Due to the high sensitivity of the fluorescent result, there is no advantage to incubating the 3M Attest Super Rapid Readout Biological Indicators 1491 and 1492V beyond 24 minutes. After activation and during incubation, the white Nonwoven Material in the biological indicator will absorb bromocresol purple indicator, the pH-sensitive indicator dye in the growth media, and appear blue. In the case of a positive control BI, a yellow color change of the growth media and/or Nonwoven Material will appear within 24 hours for 1491 BIs and within 48 hours for 1492V BIs.

Any observation of a yellow color within the vial indicates a positive result. In the case of a processed BI, a media and/or Nonwoven material color change from purple to yellow indicates a steam sterilization process failure. A negative pH color change result, i.e. media and Nonwoven Material remain purple/blue, can be assessed at 24 hours for 1491 BIs and at 48 hours for 1492V BIs.

Optional Visual pH Color Change Result - 3M Attest Rapid Readout Biological Indicator 1295

3M Attest Rapid Readout Biological Indicators 1295 for Vaporized Hydrogen Peroxide Sterilization may be further incubated for a visual pH color change. Due to the high sensitivity of the fluorescent result, there is no advantage to continuing to incubate 3M Attest Rapid Readout Biological Indicators 1295 after the fluorescent result has been determined and documented. This optional further incubation is not intended for routine use but may be used for special studies.

In the case of a positive control BI, a yellow color change of the growth media will appear within 24 hours. Any observation of a yellow color within the vial indicates a positive result. The pH color change result can be recorded and the BI discarded as soon as the yellow color is observed.

For processed 3M Attest Rapid Readout Biological Indicators 1295, the use of a humidified incubator is recommended to avoid media dry out. If a visual pH color change result is desired, transfer the processed 1295 BI to a humidified incubator after the fluorescent result has been determined in the 3M Attest Auto-reader 490. A media color change from purple to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media remains purple, can be assessed at 7 days.

FEATURES OF THE 3M ATTEST AUTO-READER 490

Audible Alarm

An audible alarm will sound, the LCD Display backlight will flash, and a (+) symbol will be displayed on the LCD panel directly under the incubating BI when a positive result is detected. An alarm will also sound if a 3M Attest Biological Indicator is removed from the 3M Attest Auto-reader 490 well before a final result has been detected. Data will be lost if the 3M Attest Biological Indicator is not replaced in the incubation well within 10 seconds.

Alarm On/Off

The audible alarm for positive results can be turned off for a particular 3M Attest Biological Indicator. This is accomplished by pressing the "Buzzer Alarm OFF" button after the alarm has sounded. This will turn off the alarm for that biological indicator only. Any positive 3M Attest Biological Indicator results subsequently detected in other wells will result in the alarm sounding. Removing the 3M Attest Biological Indicator will reactivate the alarm for that well. Pressing the "Buzzer Alarm OFF" button will not silence the alarm for error or caution code alarm warnings. To disable the alarm, see configuration section.

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Stackable

To save counter space, 3M Attest Auto-reader 490s may be stacked two high.

Information Button

When the Information button is pressed once the 3M Attest Auto-reader 490 displays two lines of text in the LCD display.

The top row indicates (left to right):

- 3M Attest Auto-reader 490 serial number
- 3M Attest Auto-reader 490 embedded software version (V X.Y.Z)
- Current date and time.

The bottom row indicates:

- The BI type configuration for wells 1-10 (i.e., Attest)

Press the information button twice and the 3M Attest Auto-reader 490 LCD panel will display:

- Top row: Dynamic IP address (if unit is connected to an Ethernet jack)
- Bottom row: Date (YYYY-MM-DD) and time that power was supplied to the 3M Attest Auto-reader

Press the information button three times and the 3M Attest Auto-reader 490 LCD display will present the unique Media Access Control (MAC) Address for the device

Temperature Monitoring

The 3M Attest Auto-reader 490 contains an internal temperature self-diagnostic program. Should the temperature fall outside the specified range (60°C +/-2°C) an error code will be displayed. The heater block temperature may be checked by depressing the Temperature Button once.

Easy Configuration

The 3M Attest Auto-reader 490 may be configured using either the buttons located on the front of the unit or using the optional 3M Attest Auto-reader 490 Web App.

Optional 3M Attest Auto-reader 490 Web App

If desired, the 3M Attest Auto-reader 490 may be connected to a PC via an Ethernet connection for additional capabilities:

- Remotely view current status of BIs in the 3M Attest Auto-reader 490 via PC.
- For full functionality, the use of an HTML5 compliant browser with JavaScript enabled is required
- Remotely view and configure 3M Attest Auto-reader 490 settings
- View the last 100 BI records
- Print out documentation using a local or networked printer and/or save documentation as an electronic file.*
- Refer to your facility's record keeping policy and procedure.

CARE AND CLEANING

CAUTION: Always unplug the 3M Attest Auto-reader 490 and allow to cool before cleaning.

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Clean the exterior surface of the 3M Attest Auto-reader 490 by wiping it with a damp cloth moistened with a mild detergent. Do not attempt to clean interior components. If additional cleaning of the wells is required, call for repair or replacement.

Note: The 3M Attest Auto-reader 490 does not contain any user serviceable parts or require any routine maintenance. Any malfunctions will be indicated by an error code. Refer to the Troubleshooting Guide section for error code explanations.

TROUBLESHOOTING GUIDE

In the event of an error code, unplug the 3M Attest Auto-reader 490 and plug it back in to initiate the self-diagnostic test. If error code re-occurs, call for repair or replacement. The LCD Display panel will indicate detected malfunctions by displaying a caution or error code under affected BI wells. The audible alarm, if active, sounds and the LCD Display backlight flashes to alert the user to all caution (except C1) and error codes. The alarm and backlight can be silenced and extinguished, respectively, by pressing the Buzzer Alarm OFF button . Note that disabling the audible alarm does not disable the flashing LCD Display backlight. Refer to the following:

Code/ Issue	Problem	Solution
C1	Unit Warm-up Caution Incubation block not at temperature.	Allow a 30 minute warm-up after 3M Attest Auto-reader 490 is plugged in. (Note: 3M Attest BI readings are disabled until incubator block reaches 60°C +/-2°C).
C2	BI Removed from Well Caution 3M Attest BI has been removed before incubation	Replace 3M Attest BI in proper well within 10 seconds to prevent loss of collected data.

Code/ Issue	Problem	Solution
	is complete.	
C4	Ambient Light Caution	Move 3M Attest Auto-reader 490 to reduced ambient light conditions.
E1	Temperature Control Error This error will occur if system is no longer able to control the temperature of the incubator block (60°C +/-2°C).	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement.
E3	Memory Error 3M Attest Auto-reader 490 unable to save/read data from memory.	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement.
E7	Communications Error This error will occur if internal communications have failed.	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement.
E8	LED Error This error occurs for a specific well if the system detects that a UV LED or Optics Detection circuit is not working.	Unplug the unit and plug it back in to clear error. If error doesn't clear, call for repair or replacement. Do not use affected incubation well.
E11	BI Inserted Incorrectly Error Incubation time remaining will not display and fluorescent readings are not taken.	Remove BI and inspect to verify cap is fully closed, glass ampoule is crushed, and growth media is in growth chamber. If necessary, repeat BI activation steps and return BI to incubation well, ensuring it is seated properly.
E12	Power Failure Error Power to Auto-reader was off for too long to provide valid result for affected BI(s).	Ensure power supply connection to Auto-reader is secure. Retest sterilizer with a new BI.
E14	Ambient Light Error (occupied well) Fluorescent result will not display.	Rely on optional visual pH color change result for BI in affected well. Move 3M Attest Auto-reader 490 to reduced ambient light conditions to avoid future occurrences. To clear E14 error code, after visual pH color change result has been documented and BI has been removed, unplug the 3M Attest Auto-reader 490H and plug it back in.
Other	Unwanted alarm sounds when positive BI detected.	Disable alarm, refer to 3M Attest Auto-reader 490 Configuration section for instructions.
Other	Alarm Sounds 3M Attest Biological Indicator has been removed from incubation/reader well prematurely.	Replace 3M Attest Biological Indicator in incubation/reader well within 10 seconds to prevent loss of data collected.
Other	LCD Panel does not display properly.	Cycle power. If problem persists call for repair or replacement.
Other	Negative Control BI BI not properly activated, Media did not contact spore carrier. Control BI was sterilized.	Ensure BI cap is closed and glass ampoule is crushed. Flick BI until media flows to spore chamber at bottom of BI. Use a non-sterilized BI as a control. Check process indicator on BI cap.
Other	Unexplained Positive BI Adhesive residue on vial. Test vial is absorbing fluorescent residue from a chemical indicator or tape.	Do not place tape or labels on vial. Place the vial so it does not come in direct contact with chemical indicators or tape.
Other	3M Attest Auto-reader 490 will not enter Configuration mode. Pushbutton Configuration may be disabled on 3M Attest Auto-reader 490 Web App.	On 3M Attest Auto-reader 490 Web App Configuration Page, go to Pushbutton Configuration and select the Enable radio button. Save changes to begin using new setting.

Code/ Issue	Problem	Solution
Other	3M Attest Auto-reader 490 LCD display flashes and asterisks (***) and/or E11 and/or C4 codes appear on display below unoccupied incubation wells.	Move 3M Attest Auto-reader 490 away from sunlight or sources of incandescent light.

OPTIONAL 3M ATTEST AUTO-READER 490 WEB APP

Web Browser Requirements

The Web Interface is designed to be viewed using an HTML5 compliant browser with JavaScript enabled (e.g., Google Chrome).

Ethernet Connection:

Connection switch/router with auto-negotiation capability to identify the Auto-reader as operating at 10 Mbps and half-duplex, or able to be manually configured as 10 Mbps and half-duplex.

Connecting to the 3M Attest Auto-reader 490 Web App

To use the 3M Attest Auto-reader 490 Web App, the 3M Attest 490 Auto-reader must be connected to the facility network using the supplied Ethernet cable. Consult with your facility IT personnel for assistance with network connectivity. Each 3M Attest Auto-reader 490 has a unique MAC Address which is presented on the unit's LCD Display when the user depresses the Info button three times. The MAC Address is also provided at the bottom of the 3M Attest Auto-reader 490 Web App Configuration page. Communication occurs using the http protocol over Port 80. Users may connect to the 3M Attest Auto-reader 490 Web App using a URL based upon either the 3M Attest Auto-reader 490 serial number or the IP Address. To view the 3M Attest Auto-reader 490 Web App from a computer on a sub-net different from that to which the 3M Attest Auto-reader 490 is connected, use the IP Address URL.

Note: To access the 3M Attest Auto-reader 490 Web App within a secure network, the PC must have access to the secure network; options include connecting via a network Ethernet cable, signing in to the facility's wireless network, or connecting via a virtual private network (VPN)..

A. Serial number URL

1. Connect the supplied Ethernet cable from the 3M Attest Auto-reader 490 to an Ethernet jack.
2. Depress the  button once to display the 3M Attest Auto-reader 490 6-digit serial number on the left side of the top row of text on the LCD display.
3. Enter `http://attestxxxxxx` where xxxxxx represents the unit's serial number (e.g. `http://attest101127`) from Step 2 in PC web browser address bar and press Enter.
4. 3M Attest Auto-reader 490 Web App Status page will display on PC. For future ease of use, add this page to your web browser "Favorites".

B. IP Address URL

The 3M Attest Auto-reader 490 Web App may also be accessed using either a dynamic or static IP address.

1. Dynamic IP Address

A dynamic IP address is assigned when the Ethernet cable is connected from the 3M Attest Auto-reader 490 to a live Ethernet jack.

Press and release the  button twice to display the 3M Attest Auto-reader 490's dynamic IP address. Enter `http://IPaddress` (e.g. `http://xxx.xxx.xxx.xxx`) in the PC's web browser address bar and press Enter. The 3M Attest Auto-reader 490 Web App Status page will display.

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
You may wish to label the Ethernet jack with the assigned address. Enter `http:// StaticIPaddress` (e.g. `http://xxx.xxx.xxx.xxx`) in the PC's web browser address bar and press Enter. The 3M Attest Auto-reader 490 Web App Status page will display.

About the 3M Attest Auto-reader 490 Web App Pages

The 3M Attest Auto-reader 490 Web App includes four pages:

a) Status page

1. User-entered nickname (optional)
2. Well Configuration
3. Well Number
4. BI Result/Remaining Incubation Time
5. Date and Time
6. 3M Attest Auto-reader 490 Serial Number

The Status page provides a real-time representation of the 3M Attest Auto-reader 490 configuration and LCD panel display. The current date and time and 3M Attest Auto-reader 490 serial number are also indicated on this page. From this page, the user may navigate to:

- the History page by clicking on History on the left of the screen.
- the Results Detail page by clicking on the well number of an occupied well (i.e. a well that currently has an incubating BI). Clicking on an unoccupied well takes the user to the History page.

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Note 1: The top four rows (highlighted in green) are populated automatically by the 3M Attest Auto-reader 490 and cannot be changed by the user. The remaining rows are completed by the user to satisfy the documentation requirements of their facility.

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The 3M Attest Auto-reader 490 stores results for the last 100 incubated BIs on embedded software. The 3M Attest Auto-reader 490 Web App History page provides a summary of these last 100 results.

d) Configuration page

This page allows the configuration of the 3M Attest Auto-reader 490's network and language settings from the 3M Attest Auto-reader 490 Web App.

Documenting BI Results using the 3M Attest Auto-reader 490 WebApp

3M									
3M™ Attest™ Auto-reader									
Status									
<u>History</u>									
ATTEST					ATTEST				
1	2	3	4	5	6	7	8	9	10
+	-	24						-	+

2. Enter the specific information for that BI (e.g. BI Type, BI Lot #, sterilizer #, load #, technician name, etc.) in the provided editable fields. Use the Notes section to capture information not prescribed in the provided row titles, for example, you may wish to use this field to record load items.

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Embedded software in the 3M Attest Auto-reader 490 saves the data from the auto-populated (i.e. green highlighted) fields. The application relies on the temporary storage in the user's web browser to save the data in the user-populated fields. In both cases, data is only stored for the last 100 records i.e., the 101st incubated BI will erase the information saved for the 1st incubated BI. To retain the Results Detail information for future reference it is best practice to print the page as soon as the BI result is available.

Using the BI History Page

To view or print an individual record, click on the ID number of the desired record to navigate to the Results Detail page for that BI

Status
History

ID	Well#	Start Time	Result
1	#3	2017-09-16 08:30:38	...
2	#9	2017-09-16 07:19:47	-
3	#2	2017-09-16 07:08:18	-
4	#1	2017-09-16 07:07:53	+
5	#10	2017-09-16 07:05:11	+
6	#9	2017-09-16 06:47:08	-
7	#10	2017-09-16 06:40:54	+
8	#2	2017-09-16 06:39:49	-
9	#1	2017-09-16 06:39:05	+
10	#10	2017-09-16 07:32:00	-
11	#9	2017-09-16 07:23:37	-
12	#10	2017-09-16 07:20:24	+

Note 1: Column 1, "ID" is the assigned ID number with ID 1 representing the newest record and ID 100 representing the oldest record. To help find a particular record, all columns are sortable, simply click once on the column heading to sort in ascending order and click a second time to reverse the order.

Note 2: The History page refreshes less frequently than the Status page. Until the BI result is available, the History page displays ... in the Results column indicating incubation is still in-process.

Using the Configuration Page

Certain set-up and configuration options are available on the password protected 3M Attest Auto-reader 490 Web App Configuration page. Ideally, the 3M Attest Auto-reader 490 should be empty (i.e. have no incubating BIs) when it is reconfigured. Navigate to the Configuration page by entering either:

1) <http://attestxxxxxx/protect/config.htm> (where xxxxxx is the 6-digit unit serial number), or

2) <http://xxx.xx.xxx.xx/protect/config.htm> (where xxx.xx.xxx.xx is the IP address)

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Language	Selected setting indicates language in which text on 3M Attest Auto-reader 490 Web App pages will appear on the computer screen.
Pushbutton Configuration	Selecting Enable ☐ allows the user to configure the 3M Attest Auto-reader 490 from the buttons on the unit itself. Selecting Disable ☐ prevents the user from being able to configure the 3M Attest Auto-reader 490 from the buttons on the unit itself.
Buzzer Configuration	Selecting Enable ☐ enables the audible alarm. The speaker icon will appear on the LCD panel. Selecting Disable ☐ disables the audible alarm. The speaker icon will disappear from the 3M Attest Auto-reader 490 LCD panel.
Attest Auto-reader 490 Nickname	Naming the 3M Attest Auto-reader 490 may be useful in facilities with multiple units. The assigned nickname will appear on the 3M Attest Auto-reader 490 Web App Status page.
Enable DHCP	Default: setting is checked indicating the Web App will use the Dynamic IP address automatically assigned when the 3M Attest Auto-reader 490 is connected to the network. Remove the check if you wish to use a Static IP address assigned by your facility IT.
IP Address	Enter Static IP address assigned by IT

Gateway	Enter Gateway
Subnet Mask	Enter Subnet Mask
Primary DNS	Enter Primary DNS
Secondary DNS	Enter Secondary DNS.
MAC Address	XX:XX:XX:XX:XX:XX (Unique MAC address assigned to device. This value cannot be changed.)
Use Computer Time	Check this box if you wish to set the clock on-board the 3M Attest Auto-reader 490 using the computer's clock.
Attest Auto-reader 490 Time	Alternatively, leave the Use Computer Time box unchecked and complete the date and time boxes to set the 3M Attest Auto-reader 490 clock.

After changing desired settings, click on Save Config button to save your changes. A new screen confirming the settings were successfully changed appears.

TROUBLESHOOTING THE 3M ATTEST AUTO-READER 490 WEB APP

Issue	Problem	Solution
Unable to access 3M Attest Auto-reader 490 Web App on networked computer	Web Pages have been disabled.	Enable Web Pages using push buttons on unit (see 3M Attest Auto-reader 490 Configuration section)
Unable to access 3M Attest Auto-reader 490 Web App on networked computer	3M Attest Auto-reader 490 and/or computer not connected to network.	Ensure 3M Attest Auto-reader 490 Ethernet cable is connected to a live jack and that connection to 3M Attest Auto-reader 490 is secure. Ensure computer is connected to network. If failure is with connection attempted using 3M Attest Auto-reader 490 serial number URL, try connecting using IP address URL.
	Auto-negotiation issue between the 3M Attest Auto-reader 490 and facility network switch.	Switch port setting must be set at 10 Mbps and half-duplex. Contact facility IT or 3M if problem persists.
Data entered by user on Results Detail page disappears	User entered data is written to temporary web browser storage. Browser may have cleared data.	Set web browser to accept (allow or enable) cookies. As a best practice, print or save Results Detail record as soon as BI result is final.
BI result displayed on 3M Attest Auto-reader 490 LCD panel does not appear on Results Detail page	The Results Detail page has not refreshed yet.	Use your browser to refresh the Results Detail page, or wait up to 60 seconds for the page to automatically refresh.

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Maintenance and repair

Repairs must be done by a qualified service person. Attest incubators that require service or repair should be returned, contact the local subsidiary. The 3M Attest Incubator is serviceable electronics that does not have a designated shelf life due to the ability to repair and replace parts. The duration of the warranty period is 12 months.

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The screenshot shows the 3M Attest Auto-reader interface. At the top left is the 3M logo. Below it is a black bar with the text '3M™ Attest™ Auto-reader'. On the left side, there are two links: 'Status' and 'History'. The main display area shows two 'ATTEST' test results. The first result is a 2x10 grid of numbers and symbols. The second result is a 2x10 grid of numbers and symbols. The first result has a circled '2' in the second position of the first row. An arrow points from the bottom left towards the circled '2'. Below the grids, there is a clock icon and the time '2017-09-16 08:31', and a serial number 'SN 100157'. At the bottom, there is a copyright notice: 'Copyright © 2017 3M, Inc.'

ATTEST					ATTEST				
1	2	3	4	5	6	7	8	9	10
+	-	24						-	+

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