This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

Apostille

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1. District of Columbia, United States of America
KAMAL S. AHMAD
2. This public document has been signed by _____
NOTARY PUBLIC IN AND FOR THE DISTRICT OF COLUMBIA
3. acting in the capacity of _____
KAMAL S. AHMAD, NOTARY PUBLIC IN AND FOR THE
DISTRICT OF COLUMBIA
4. bears the seal/stamp of _____

CERTIFIED

5. at Washington, D.C.

6. the 22, day of JULY 2020

7. by Secretary of the District of Columbia

8. No. 573472

9. Seal/Stamp

10. Signature:

Kimberly A. Bassett

Kimberly A. Bassett
Secretary of the District of Columbia

Self-Truthfulness Declaration

Manufacturer


Siemens Medical Solutions USA, Inc.
22010 S.E. 51st Street
Issaquah, WA 98029, USA

As a responsible individual authorized to represent and act on behalf of Siemens Medical Solutions USA, Inc. Ultrasound Business Unit, I hereby certify to the best of my knowledge that the documents provided as below for the following product used in Russia is genuine, complete and valid:

User Manual Supplement for Russia

ACUSON Sequoia Diagnostic Ultrasound System with Accessories

Signed in July 10, 2020


HyunJung Lee

Regulatory Technical Specialist
Siemens Medical Solutions USA, Inc.

**I Certify That This Is An
Original Document**

District of Columbia

Subscribed and sworn to before me,

this 21st day of July, 2020

Kamal S. Ahmad, Notary Public, D.C.

My Commission expires November 14, 2024

Siemens Medical Solutions USA, Inc.
22010 S.E. 51st Street, Issaquah, WA 98029, USA

User manual supplement for Russia

ACUSON Sequoia Diagnostic Ultrasound System with Accessories

1. The name of the Medical Device

ACUSON Sequoia Diagnostic Ultrasound System with Accessories

2. Information about Manufacturer of the Medical Device

2.1. Designated name of the manufacturer (full name, and short name if applicable), address of the manufacturing site, post address (if differs from the official company registration):

Legal Manufacturer

Siemens Medical Solutions USA, Inc.

22010 S.E. 51st Street, Issaquah, WA 98029, USA

Phone: +1-888-826-9702

www.siemens-healthineers.com

2.2. Addresses of the additional manufacturing sites (if there are any)

1) Jabil, Inc. (3800 Giddings Road, Auburn Hills, Michigan, 48326, USA)

2) Siemens Medical Solutions USA, Inc. (22010 S.E. 51st Street, Issaquah, WA 98029, USA)

2.3. The name and address of the local representative in Russia

Siemens Healthcare LLC

Dubininskaja street 96,

115093 Moscow, Russia

Tel.: +7 495 737 1252

www.siemens-healthineers.com/ru/

3. Designation of the medical device: Purposes for use (Medical Device operation data)

3.1. Intended use

The ultrasound imaging system is intended for the following applications: Fetal, Abdominal, Pediatric, Small Parts, OB/GYN (useful for visualization of the ovaries, follicles, uterus and other pelvic structures), Cardiac, Pelvic, Vascular, Musculoskeletal, and Peripheral Vascular applications.

The system also provides the ability to measure anatomical structures fetal, abdominal, pediatric, small organ, cardiac, transrectal, transvaginal, peripheral vessel, musculoskeletal, and calculation packages that provide information to the clinician that may be used adjunctively with other medical data obtained by a physician for clinical diagnosis purposes.

One of function of using the ultrasound system is biopsy control for determination of primary etiology of a disease and extent of defeat. The biopsy procedures is possible in all types of ultrasound examinations for more exact intake of material.

3.2. Indications

The ultrasound imaging system is intended for the following applications: Fetal, Abdominal, Pediatric, Small Parts, OB/GYN (useful for visualization of the ovaries, follicles, uterus and other pelvic structures), Cardiac, Pelvic, Vascular, Musculoskeletal, and Peripheral Vascular applications.

The system also provides the ability to measure anatomical structures fetal, abdominal, pediatric, small organ, cardiac, transrectal, transvaginal, peripheral vessel, musculoskeletal, and calculation packages that provide information to the clinician that may be used adjunctively with other medical data obtained by a physician for clinical diagnosis purposes.

One of function of using the ultrasound system is biopsy control for determination of primary etiology of a disease and extent of defeat. The biopsy procedures is possible in all types of ultrasound examinations for more exact intake of material.

3.3. Possible side effects

WARNING: Siemens makes every effort to manufacture safe and effective transducers. You must take all necessary precautions to eliminate the possibility of exposing patients, operators, or third parties to hazardous or infectious materials. These precautions should be considered in the use of any application that may indicate the need for such care, and during endocavity or intraoperative scanning; during biopsy or puncture procedures; or when scanning patients with open wounds.

WARNING: To eliminate the possibility of exposing patients, operators, or third parties to hazardous or infectious materials, always dispose of hazardous or infectious materials according to local, state, and regional laws and regulations.

WARNING: There have been reports of severe allergic reactions to medical devices containing latex (natural rubber). Health care professionals are advised to identify latex-sensitive patients and to be prepared to treat allergic reactions promptly. For additional information in the U.S.A., refer to FDA Medical Alert MDA91-1.

WARNING: Ultrasound energy is transmitted more effectively through water than through tissue. When using a standoff device of any kind, for example, a gel pad, the actual mechanical and thermal indices, MI and/or TI, may be higher than indicated in the output display on the system.

The assessment of the biological effects of diagnostic ultrasound on humans is a subject of ongoing scientific research. This system, and all diagnostic ultrasound procedures, should be used for valid reasons, for the shortest possible period of time, and at the lowest mechanical and thermal indices necessary to produce clinically acceptable images.

According to the ALARA (As Low As Reasonably Achievable) principles, acoustic output should be set to the *lowest level required to satisfactorily perform the examination*.

The ultrasound imaging system complies with the standards of the American Institute of Ultrasound in Medicine (AIUM), the National Electrical Manufacturer's Association (NEMA), the guidelines of the United States Food and Drug Administration (FDA), and the guidelines of the International Electrotechnical Commission (IEC) in terms of safety and acoustic output levels. The ultrasound output levels are stated to permit the user to critically evaluate the system settings in the event of new research findings being announced.

Note on Fetal Examinations

The following recommendation is excerpted from the National Institute of Health in the United States of America. Consensus Statement on the Use of Ultrasound Imaging During Pregnancy, Volume 5, No. 1, based on the recommendation issued at the Health Consensus Development Conference, February, 1984:

Ultrasound examination in pregnancy should be performed for a specific medical indication. The data on clinical efficacy and safety do not allow a recommendation for routine scanning at this time.

Ultrasound examination performed solely to satisfy the family's desire to know the fetal sex, to view the fetus, or to obtain a picture of the fetus should be discouraged. In addition, visualization of the fetus solely for educational or commercial demonstrations without medical benefit should not be performed.

In August 1994, the Food and Drug Administration (FDA) notified the medical community and the ultrasound industry regarding its concerns about the misuse of diagnostic ultrasound equipment for non-medical purposes, and to discourage patients from having sonograms for non-medical reasons.

The American Institute of Ultrasound in Medicine (AIUM) has also advocated the responsible use of diagnostic ultrasound for all fetal imaging (August 2005).

3.7. Additional information about Conditions for use (in hospitals / in fields / at home etc.)

Transducers

	During Operation
Atmospheric pressure:	700 hPa to 1060 hPa
Relative humidity:	10% to 80%, non-condensing
Temperature for 9C3, 5C1, DAX, 18L6, 14L5, 10L4, 8V3, 5V1, 4V1, CW2:	+10°C to +40°C
Temperature for 9EC4 transducer:	+10°C to +42°C

Needle guides for 9C3, 5C1, DAX, 18L6, 14L5, 10L4, 4V1 transducers

	During Operation
Temperature:	Room Temperature No specific requirements
Atmospheric pressure:	No specific requirements
Relative humidity:	No specific requirements

Disposable needle guides for 9EC4 transducers

	During Operation
Temperature:	Room Temperature No specific requirements
Atmospheric pressure:	No specific requirements
Relative humidity:	No specific requirements

Reusable needle guides for 9EC4 transducers

	During Operation
Temperature:	Room Temperature No specific requirements
Atmospheric pressure:	No specific requirements
Relative humidity:	No specific requirements

4. Classification of the Medical Device:

4.1 Type of the Medical Device in accordance with the Russian Nomenclature Classification

260250

4.2. Class of the potential risk of the medical device application

Class 2a

4.3 Class of the Medical Device in accordance with IEC 60601-1

Classifications	Contents
Type of protection against electrical shock	Class I, external powered
Degree of protection against electrical shock	– Type BF applied part for endocavity, linear, curved, and phased array transducers – Type B applied part for Fusion transducer sensors and needle tracking sensors – Type BF defibrillation-proof applied part for ECG
Degree of protection against harmful ingress of water	IPX8 for transducers, footswitch IPX0 - system
Degree of safety of application in the presence of a flammable anesthetic material with air or with oxygen or nitrous oxide	Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
Mode of operation	Continuous operation

4.4 Software Risk classification / Level of concern

Risk classification in accordance with IEC 62304: B

5. Technical Specifications of the Medical Device

5.1. Photos of the device

Photos of the system:

5.2. The list of components of the medical device

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15. Инструкция по эксплуатации - 1 шт. (но не более 5 шт. при необходимости)	User documentation – 1 pcs.(no more than 5 pcs. if necessary)	• Conventions and typographical conventions used in the manuals • Intended audience • Introduction to the ultrasound system • Safety and care information for the ultrasound system and compatible transducers • Procedures for setting up and preparing the system for use • Procedures for registering a patient and activating an operating mode • Procedures for acquiring, optimizing, annotating, measuring, printing, and recording images • Technical description of the ultrasound system • Acoustic output data
16. Приложение к инструкции по эксплуатации для России – 1 шт.	Addendum to the user manual for Russia – 1 pcs.	Includes additional information about system
17. Руководство по расширенной визуализации – 1 шт.	Advanced Imaging Manual – 1 pcs.	This manual is a companion to the Instructions for Use publication and contains additional instructions for use for the safe and proper use of advanced imaging features on the ultrasound system. • Description of the biopsy function • Procedures for specialty imaging • Panoramic imaging • Contrast agent imaging • Comparing images from multiple modalities • Virtual Touch applications • Procedures for reviewing images and clips • Procedures for the physiologic function • Fusion imaging • Velocity vector imaging
18. Дополнение к руководству пользователя и справочного руководства – 1 шт. (но не более 5 шт., при необходимости).	Supplement to the User and Reference Manual – 1 pcs. (no more than 5 pcs.).	The information in this supplement provides updates to your existing user manuals. For a complete set of operating instructions, you must use this supplement and your ultrasound product's user and reference manuals.
Принадлежности:	Accessories:	
1. Буквенно-цифровая клавиатура	Alphanumeric Keyboard	Language-specific alphanumeric keyboard available in: English, German, French, Spanish, Italian, Danish, Finnish and Swedish, Norwegian

2. Лицензионный ключ для функции визуализации с контрастным агентом	License key for Contrast Agent Imaging Option	Supports the following contrast agent destruction techniques: Burst, Flash Sequencing - Compatible transducers: 10L4 Supported studies: Abdomen, Arterial, Breast, Carotid, Musculoskeletal, Testis, Thyroid, Venous, General - Compatible transducers: 14L5 Supported studies: Breast, Carotid, Musculoskeletal, Testis, Thyroid, Venous, General - Compatible transducers: 18L6 Supported studies: Breast, Musculoskeletal, Testis, Thyroid, General - Compatible transducers: 5C1, DAX, 9C3, 4V1 Supported studies: Abdomen, Gynecology, General - Compatible transducers: 9EC4 Supported studies: Gynecology, Prostate, General - VueBox Option - Requires clips acquired during contrast agent imaging and stored in DICOM format - Enables compatibility with the VueBox software application from Bracco. The software is used to perform qualitative and quantitative analysis of linear data. For information regarding the software application, refer to the manufacturer's instructions. - Compatible transducers: 10L4, 14L5, 18L6, DAX, 5C1, 9C3, 4V1, 9EC4
3. Лицензионный ключ для функции визуализации деформации Virtual Touch.	License key for Virtual Touch Strain Imaging Option	- Provides a qualitative representation of relative tissue stiffness for the region of interest - Strain Ratio provides a quantifiable method to compare the relative stiffness of tissue within two user-selectable regions of interest - Compatible transducers: 9EC4 Supported studies: Prostate, General - Compatible transducers: 10L4, 14L5, 18L6 Supported studies: Breast, Thyroid, General
4. Лицензионный ключ для функции точечной эластографии сдвиговых волн Virtual Touch.	License key for Virtual Touch Point Shear Wave Option	- Measures tissue shear velocity (Vs) and elasticity (E) for a selected region of interest using acoustic radiation force impulse (ARFI) technology of diagnostic ultrasound to induce tissue displacement. - Provides measurement labels for sites, lesions, and liver segments - Compatible transducers: 4V1, 5C1, DAX

		Supported studies: Abdomen, General Compatible transducers: 10L4 Supported studies: Abdomen
5. Лицензионный ключ для функции эластографии сдвиговых волн Virtual Touch.	License key for Virtual Touch Shear Wave Option	- Uses acoustic radiation force impulse (ARFI) technology of diagnostic ultrasound to induce tissue displacement - Qualitatively depicts shear velocity for a selected region of interest on a 2D image - Provides quantitative measurements of shear velocity (Vs) and elasticity (E) for selected points within the region of interest - Compatible transducers: 5C1, DAX Supported studies: Abdomen, General - Compatible transducers: 10L4 Supported studies: Breast, Thyroid, General
6. Лицензионный ключ для функции визуализации в режиме слияния (Fusion).	License key for Fusion Imaging Option	- Fusion imaging aligns reference data with a real-time ultrasound image for diagnosis and interventional procedures - Supported reference data: computed tomography, magnetic resonance imaging - Compatible transducers: 4V1, 5C1, DAX - Supported studies: Abdomen, General - Compatible transducers: 10L4 Supported studies: Musculoskeletal, General - Fusion supports the following components – Ascension Technologies Corporation driveBAY 2 electromagnetic tracking system – Ascension Technologies Corporation Mid-range Transmitter (for the driveBAY 2) – GCX pole to support mid-range transmitter – CIVCO General Purpose 8 mm sensors Note: For information on the installation, care, and operation of hardware components, refer to the manufacturer's operating instructions that accompanied the device.

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5.3. Additional Technical description of the main system's components

Tolerances: $\pm 5\%$ for all the technical specifications if not indicated another information.

Contents	Specification
Maximum MD movement force	< 150N

System software

Contents	Specification
Software version	VA10
Platform	Windows
Software release date	May 3, 2018

Wireless Network Connections

Contents	Specification
Network standard:	802.11a 802.11b 802.11ac 802.11g 802.11n
Frequency bandwidth:	2.4 GHz 5 GHz
Authentication:	WPA WPA2 WPA PSK WPA2 PSK Open
Encryption:	None TKIP AES
Extensible Authentication Protocol (EAP):	PEAPv0 (PEAP-MSHCAPv2) TLS

Power cord

Contents	Specification
Length	4,5 m

Transducers specifications:

Transducer	Dimensions (mm)	Weight (g)	Number of elements	Footprint (mm)	Maximum Field of View, ° or mm	Cable length, m	Radius of curvature (mm)
10L4	138,6 x 45,2 x 14,7	723.2	192 x 2	49,25 x 18,85	38,2 mm	≤ 2.3 m	NA
14L5	138,6 x 45,1 x 8,3	726.9	192 x 2	49,58 x 12,89	38,2 mm	≤ 2.3 m	NA
18L6	112,0 x 62,0 9,2	761.8	576	69,22 x 16,48	57,5 mm	≤ 2.3 m	NA
DAX	140,23 x 54 x 27,2	878.0	96 x 2	57,74 x 30,23	50°	≤ 2.9 m	56,7
5C1	141,37 x 58,79 x 13,8	756.9	180	63,3 x 18,2	72°	≤ 2.3 m	46,0
9C3	128,7 x 63,3 x 14	780.4	384	69,56 x 20,47	78,6°	≤ 2.3 m	45

4V1	139 x 31 x 16	639.1	128	35,5 x 20,2	90°	≤ 2.3 m	NA
5V1	124,8 x 23,4 x 15,2	622.6	96	27,2 x 18,7	90°	≤ 2.2 m	NA
8V3	124,2 x 22,5 x 11,5	643.6	128	26,9 x 16,6	90°	≤ 2.3 m	NA
9EC4	340 x 20 x 21,7	700	192	17,0 x 22,0	176°	≤ 2.2 m	10,8
CW2	20 x 65 x 104.14 (17 – diameter)	100	2	NA	NA	≤ 2.3 m	NA

Technical specifications of the components

Needle guide for 9C3 transducer

Contents	Specification
Dimensions	1.10 x 3.36 x 2.11 inch
Weight	17g
Производитель	CIVCO Medical Solutions

Needle guide for 5C1 transducer

Contents	Specification
Dimensions	1.82 x 3.11 x 1.09 inch
Weight	9,4 g
Manufacturer:	CIVCO Medical Solutions

Needle guide for DAX transducer

Contents	Specification
Dimensions	2.11 x 3.16 x 1.09 inch
Weight	10.7 g
Manufacturer:	CIVCO Medical Solutions

Needle guide for 18L6 transducer

Contents	Specification
Dimensions	0.865 x 4.028 x 1.913 inch
Weight	15g
Manufacturer:	CIVCO Medical Solutions

Needle guide for 14L5 transducer

Contents	Specification
----------	---------------

Dimensions	33.27 x 64.26 x 25.65 mm
Weight	6,5 g
Manufacturer:	CIVCO Medical Solutions

Needle guide for 10L4 transducer

Contents	Specification
Dimensions	43.42 x 63.25 x 26.42 mm
Weight	8.2g
Manufacturer:	CIVCO Medical Solutions

Reusable Needle guide for 9EC4 transducer

Contents	Specification
Dimensions	29.72 x 156.97 x 27.18 mm
Weight	55g
Manufacturer:	CIVCO Medical Solutions

Disposable Needle guide for 9EC4 transducer

Contents	Specification
Dimensions	22.35 x 156.46 x 26.67 mm
Weight	13.7g
Manufacturer:	CIVCO Medical Solutions
Shelf life	3 years
Sterilization type	EO
Disposability	Disposable

Needle guide for 4V1 transducer

Contents	Specification
-----------------	----------------------

Dimensions	32.00 x 65.33 x 47.57 mm
Weight	18g
Manufacturer:	CIVCO Medical Solutions

BW printer

Recommended BW printer - UP-X898MD, manufactured by Sony

Contents	Specification
Dimensions	154 × 88 × 240 mm
Weight	2,5 kg
Power voltage	AC 100V to 240V, 50/60 Hz
Input current	1,3 A to 0,6 A
Resolution	325 dpi
Printing speed	Approx 1,9 sec / image at standard settings (High speed print mode) Approx 3,3 sec / image at standard settings (Normal speed print mode)
Gradation	256-level gradation (8-bit)
Memory	4096 × 1280 dots (for digital)

Color printer

Recommended Color printer - UP-D25MD, manufactured by Sony

Contents	Specification
Dimensions	212 × 98 × 398 mm
Weight	5,5 kg
Power requirements	AC 100 to 240V, 50/60 Hz
Input current	1,7 – 1,0 A
Resolution	423 dpi
Printing speed	21L – approx. 29 sec; 21S – approx. 19 sec; 24LA – approx. 36 sec (tentative); 24SA – approx. 25 sec (tentative)
Gradation	256-level gradation (8-bit)

Configurable three-pedal footswitch

Contents	Specification
Manufacturer	Kinesis
Dimensions, mm	27.3 x 17.8 x 3.2 cm
Weight, g	1.4 kg
IPXX	IPX8
Power supply voltage	+5V +/- 0.5V, 50mA max.
Withstand voltage	AC 500V 60sec.
Cable length	~8ft (~2.4m)

Accessory shelf

Content	Specification
Dimensions	274(W) x 235.72(L) x 65.23(H)mm
Weight	0,5 kg
Maximum load	5 kg

Universal Video Converter

Contents	Specification
Manufacturer	Kortek Corporation
Brand name	UVC-300K
Rated power	12V DC, Max 100 mA
Power consumption	Max 12W
Video Input (Resolution)	1920 x 1080 HD pixel format
Video Output (Resolution)	DVI/ D-Sub : 1024 x 768 60Hz NTSC 720 x 480 at 15.7kHz/60Hz PAL 720 x 575 at 15.66kHz/50Hz - One of 2 S-Video types selected by user using DIP switch selection.
Weight	Not to exceed 0,9 kg based on the main module
Dimensions	210 (W) x 50 (H) x 180 (D) mm

ECG cable

Manufacturer:	Shenzhen Med-Link Electronics Tech Co., Ltd
Type of connectors	EU (IEC)

Number of leads	3
Cable length	3050±50mm
Leads length	1000±50mm

eSle link accessories:

Contents	Specification
Camera:	
Model:	Logitech C930E
Manufacturer	Logitech
Video performance	operation with 1080P video at 30 frames-per-second
USB headset	
Model:	JABRA EVOLVE 30 UC MONO
Manufacturer	Jabra
Performance	USB 2.0 compliant

Accessories for the visualization function in merge mode (Fusion):

- Electromagnetic tracking system;

Contents	Specification
Manufacturer	Ascension
Dimensions, mm	180 mm x 147 mm x 41 mm
Weight, kg	0.84 kg
Power	Molex Power Connector; +12 V: 1.6 A nominal, 2.9 A maximum; +5V: 600 mA nominal

- Transmitter SV - band;

Contents	Specification
Manufacturer	Ascension
Dimensions, mm	96 mm x 96 mm x 96 mm
Weight, kg	2.3 kg
Minimum operating range:	5.7 cm from any face of the MRT transmitter housing

Accuracy	less than 1.4 mm RMS, 0.5 degrees RMS
Measurement/Update rate	The transmitter will operate at system measurement rates between 20 and 250 Hz.

- Stand GCX to install the transmitter CB - band;

Contents	Specification
Manufacturer	GCX Mounting Solutions
Height	Min. 122 cm
Load	Min. 9.5 kg

- General Purpose Sensor;

Contents	Specification
Manufacturer	CIVCO Medical Solutions
Dimensions, mm	Sensor Housing: 12.2 x 12.4 x 49.5 (Cable Length 17.78 cm)
Weight, kg	64.76 g

Additional technical specifications:

Technical specification	Value
Maximum allowable time for setting the operating mode	160 seconds as booting time
Corrected level of noise (sound power)	NA, the patient and user is not exposed to acoustic noise
Effective radiated power of a wireless transmitter of the device	The SparkLan wireless transmitter is limited to 15dBm with a tolerance of +/- 2dBm

5.4. Electromagnetic Compatibility

The electromagnetic compatibility (EMC) testing of this system was performed according to the international standard for EMC of medical devices (IEC 60601-1-2). This International Electrotechnical Commission (IEC) standard has been adopted in Europe as the European Norm (EN) standard (EN 60601-1-2). The ultrasound system from Siemens meets the requirements of the standard. Compliance with test limits does not guarantee that a particular installation will be free from electromagnetic interference.

Avoiding Electromagnetic Interference

WARNING: Operating the ultrasound system in close proximity to other equipment can cause reciprocal interference. You should observe and ensure normal operation of the ultrasound system and other equipment.

Medical devices either generate (emit) or receive electromagnetic interference. The EMC standards describe tests for both emitted and received interference. Emission tests address interference generated by the device being tested. The ultrasound system from Siemens meets the standard for radiated emissions.

Radio frequency (RF) is a form of electromagnetic energy tested by the EMC standards. Ultrasound systems and ultrasound transducers operate in the range of RF frequencies and are therefore susceptible to interference generated by other RF energy sources. Sources of RF interference include medical devices, information technology products, and radio/television transmission towers. To reduce interference from an RF energy source:

- Increase the separation (distance) between the ultrasound equipment and the interfering device (RF energy source).
- Connect the ultrasound equipment and the interfering device (RF energy source) to separate outlets on separate circuits, if possible.

In some environments, it is difficult to identify the source of radiated interference. To identify the source of the interference, determine the answers to these questions and then contact your service representative with the information you have collected.

- Is the interference intermittent or constant?
- Is the interference present with only one transducer or with several transducers?
- Do two different transducers operating at the same frequency have the same problem?
- Is the interference present if the ultrasound system is moved to a different location in the facility?
- Is the interference mode-dependent? That is, does the interference occur only in one mode, such as during Doppler?
- Can the EMC coupling path be attenuated? For example, does placing a transducer close to an ECG cable increase electromagnetic interference? Does moving the cable or other medical equipment away from the transducer result in reduced electromagnetic interference?

Electrosurgical Units

Electrosurgical units (ESUs) and other devices intentionally introduce radio frequency electromagnetic fields or currents into patients. Because ultrasound imaging frequencies are coincidentally in this radio frequency range, ultrasound transducer circuits are susceptible to radio frequency interference. When an ESU is in use, expect the ultrasound image to experience severe interference that may make the image unusable. This interference stops as soon as the ESU is powered off.

Restrictions for Use

If RF interference is detected, then the physician must determine if an artifact caused by the interference will negatively impact image quality and the subsequent diagnosis.

Immunity Level Test Results

Immunity is defined in the standard as the ability of a system to perform without degradation in the presence of electromagnetic disturbance. The EMC standards require manufacturers of patient-coupled equipment to specify immunity levels for their systems. The standards recognize

that ultrasound equipment is designed to receive and amplify low-level signals in the same bandwidth as the interference. Therefore, it is reasonable to expect image noise during electromagnetic disturbance.

EMC Note: Operating the ultrasound system in close proximity to sources of strong electromagnetic fields, such as radio transmitter stations or similar installations, as well as portable and mobile RF (radiated frequency) may lead to interference visible on the monitor screen. However, the device has been designed and tested to withstand such interference and will not be permanently damaged.

When comparing immunity levels of different ultrasound systems, recognize that although the EMC standard does define the test methodology, it does not specify the criteria to assess degradation. Degradation assessment may vary by manufacturer.

A qualitative assessment of degradation in image quality is subjective. Noting when the first sign of an artifact is seen in the image minimizes the issue of subjectivity and also provides for stringent test results.

The results of emissions testing and immunity testing are provided in the accompanying tables. The standards used in testing are also provided in the tables. Testing was performed on a typical ultrasound system configuration. The values in the compliance level column of the tables specific to transducers is the level at which the RF interference is no longer discernable from the ambient background noise.

Transducers, Accessories, and Cables

WARNING: Accessory equipment connected to the analog and digital interfaces must be certified according to the respective EN or IEC standards (for example, EN 60950 or IEC 60950 for data processing equipment and EN 60601-1 or IEC 60601-1 for medical equipment). Anyone who connects additional equipment to any of the signal input or signal output ports configures a medical system and is therefore responsible that the system complies with the requirements of the system standards EN 60601-1 or IEC 60601-1. Siemens can only guarantee the performance and safety of the devices listed in the Instructions for Use. If in doubt, consult the Siemens service department or your local Siemens representative.

WARNING: You must only use the transducers, accessories, cables, and replacement parts for internal components specified by Siemens to reduce the risk of increased RF (radio frequency) emissions or decreased immunity of the ultrasound system.

Guidance and manufacturer's declaration — electromagnetic emissions for the ACUSON Sequoia system

The ACUSON Sequoia system is intended for use in the electromagnetic environment specified below. The customer or the user of the ACUSON Sequoia system should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The ACUSON Sequoia system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any

		interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The ACUSON Sequoia system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration — electromagnetic immunity for the ACUSON Sequoia system independent of transducer type

The ACUSON Sequoia system is intended for use in the electromagnetic environment specified below. The customer or the user of the ACUSON Sequoia system should assure that it is used in such an environment.			
Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle	System may shut down and the mains circuit breaker could be switched to the off position (⏻ = OFF).	Ensure the mains circuit breaker is in the on position (I = ON) and then use the power on sequence to reboot the system. Mains power quality should be that of a typical commercial or hospital environment. If the user of the ACUSON Sequoia system requires continued operation during power mains interruptions, it is recommended that the

			ACUSON Sequoia system be powered from an uninterruptible power supply or a battery.
	40% U_T (60% dip in U_T) for 5 cycles	40% U_T (60% dip in U_T) for 5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ACUSON Sequoia system requires continued operation during power mains interruptions, it is recommended that the ACUSON Sequoia system be powered from an uninterruptible power supply or a battery.
	70% U_T (30% dip in U_T) for 25 cycles	70% U_T (30% dip in U_T) for 25 cycles	
	<5% U_T (>95% dip in U_T) for 5 seconds	System goes to normal power down condition.	For extended mains power outages without the use of an uninterruptible power supply or a battery, the system will shut down normally. Use the power on sequence to reboot the system.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration — electromagnetic immunity for the ACUSON Sequoia system

The ACUSON Sequoia system is intended for use in the electromagnetic environment specified below. The customer or the user of the ACUSON Sequoia system should assure that the system is used in such an environment.			
Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the ACUSON Sequoia system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			Recommended separation distance (<i>d</i>) $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ACUSON Sequoia system is used exceeds the applicable RF compliance level above, the ACUSON Sequoia system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as changing transducers or transducer operating frequency, or re-orienting or relocating the ACUSON Sequoia system. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than the compliance level 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the ACUSON Sequoia system

The ACUSON Sequoia system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ACUSON Sequoia system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment

(transmitters) and the ACUSON Sequoia system as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W (watts)	Separation distance according to frequency of transmitter m (meters)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

6. Sterility (if applicable).

Diagnostic ultrasound system ACUSON Sequoia – non-sterile.

Needle guide for 9C3 transducer; needle guide for 5C1 transducer; needle guide for DAX transducer; needle guide for 18L6 transducer; needle guide for 14L5 transducer; needle guide for 10L4 transducer; needle guide for 4V1 transducer – non sterile, must be sterilized before using according to instructions for use.

Reusable needle guide for 9EC4 transducer – non sterile, must be sterilized before using according to instructions for use.

Disposable needle guide for 9EC4 transducer – sterile (EO).

7. Additional Transportation requirements:

Needle guides for 9C3, 5C1, DAX, 18L6, 14L5, 10L4, 4V1 transducers

	During Storage or Transportation
Temperature:	+11°C to +29°C
Relative humidity:	No specific requirements
Atmospheric pressure:	No specific requirements

Disposable needle guides for 9EC4 transducers

	During Storage or Transportation
Temperature:	+11°C to +29°C
Relative humidity:	No specific requirements
Atmospheric pressure:	No specific requirements

Reusable needle guides for 9EC4 transducers

	During Storage or Transportation
Temperature:	+11°C to +29°C
Relative humidity:	No specific requirements
Atmospheric pressure:	No specific requirements


8. Labeling.

8.1 Symbols

Symbol	Explanation
	AC (alternating current) voltage source
	Caution: Risk of electric shock (black and yellow illustration)
	Defibrillation-proof type BF applied part
	Type BF applied part
	Consult Instructions for Use
	It is mandatory you refer to the manual. (blue illustration)
	Caution (black illustration)
	General Warning (yellow and black illustration)
	Attention (black illustration)
	Partial power on/off control
	ON position of the mains switch (circuit breaker on the ultrasound system)

Symbol	Explanation
	OFF position of the mains switch (circuit breaker on the ultrasound system)
	Equipotential terminal.
	Protective Earth Ground
	Signal input
	Signal output
	USB Connection
	Direct current (DC) power input connection
	Indicator light
	Ethernet Connection
	Auxiliary continuous wave Doppler transducer port
Aux 1	Auxiliary in/out port
ECG	Electrocardiogram (EKG)
	ECG signal input connection
	Auxiliary in/out port
	Unlocked position
	Locked position
	Swivel
	Left-right rotation adjustment
	Height adjustment
	Swivel lock engaged
	Lock engaged
	Wheel lock engaged

Symbol	Explanation
	Wheel lock released
	Use only in dry conditions Do not install a wet air filter.
	Location of air filter
PUSH	Releases or engages the locking mechanism when pushed
IPX8	Protected against the effects of continuous immersion in water
	Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of, a physician.
	Safe working load "XXX" represents the weight of the system unpacked and ready for operation with documentation devices installed Additional objects or devices on the system cannot exceed the safe working load
	The product must be properly disposed of in accordance with local, state, and regional laws and regulations. Products bearing this symbol are subject to the European Community directive 2002/96/EC on waste electrical and electronic equipment (WEEE), amended by directive 2003/108/EC. For collection and disposal of the product, its components, or its accessories, contact your local Siemens representative.
	Environmentally friendly use period
	Do not exert excessive force when moving or positioning the ultrasound system. (Black and red illustration)
	Do not lean. (red and black illustration)
	Shelf or compartment weight restriction "XX" represents the maximum weight restriction
	Pinch hazard
	Manufacturer's declaration of product compliance with applicable EEC directive(s) and the European Notified Body
	Manufacturer's declaration of product compliance with applicable EEC directive
	UL symbol for listing as recognized components for Canada and United States of America


Symbol	Explanation
	UL classified symbol for Canada and United States of America
	Authorized representative in the European community
	EAC (Eurasian Conformity mark) The Eurasian Conformity mark is a product marking to indicate that the product conforms to all Technical Regulations of the Customs Union assessment procedures
	Equipment complies with the standard and specifications published by IMDA (Info-communications Media Development Authority of Singapore) "XXXXXX" represents the dealer's license number
	TÜV Rheinland INMETRO Certification Mark (Brazil)

Symbol	Explanation
 N232	Product complies with applicable Australian or New Zealand standards for EMC and radio-communications.
	DEMKO-Denmark approval mark
	National certification for Korea
	Ukrainian conformity mark
	The monitor is compliant with the limits for a Class B digital device, pursuant to Part 15 of the Federal Communication Commission (FCC) rules.
	Location of device manufacturing: United States of America
	Location of device manufacturing: Republic of Korea
	Location of device manufacturing: Japan
IVK	(Installierte Volumen Komponente) Identifier of selected system components or parts for product traceability
System IVK	(System Installierte Volumen Komponente) Identifies the system for product traceability
	Barcode with product identification information according to the following barcode specification: GS1 DataMatrix
	Manufacturer
 20XX	Date of Manufacture symbol with the date below
GTIN	Global Trade Item Number
Rev REV	Model revision
Serial SERIAL	Serial number
Count COUNT	Quantity of product
Model MODEL	Siemens model number
REF	Device part number
Origin US ORIGIN US	Manufacturer country of origin: United States of America
Origin KR ORIGIN KR	Manufacturer country of origin: Republic of Korea


Symbol	Explanation
Origin CN ORIGIN CN	Manufacturer country of origin: China
Origin JP ORIGIN JP	Manufacturer country of origin: Japan
	Storage temperature range (example)
	Weight and size "XXX" represents the weight of the system packaged for shipping "XX" represents the size of the box in inches
	Do not stack
	Stacking limitation "n" indicates the number of items that can be stacked
	Recycle
	The wood packaging material for shipping the ultrasound system complies with the International Standards for Phytosanitary Measures No. 15. (example)
	Do not expose to direct sunlight
	Storage atmospheric pressure range (example)
	Storage humidity range (example)
	Indicates this side up
	Do not allow to get wet
	Fragile. Handle with care.
	Refer to the operator's manual for information about compatible needle guides.
	Latex
	Don't use in case of packing damage
	EO sterilization
	Non sterile

Symbol	Explanation
	Single use
	Device class II

8.2 Labels positions


Example of the location of the system label.

Example of ultrasound system label with certification labels.


- 1 Identification label
- 2 Country-specific labels

Example of identification label.

- 1 Product name
- 2 System input power requirement



Additional Label location:


1 – Label pinch weight

2 – Label pinch



3 – Label system rear

4 – Label weight limitation


5 – Label weight white

Example of labels on a transducer.

Real photo of the main system's label:

Artwork of the medical device labeling in Russian language:

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Note: The legal manufacturer address has been changed to "22010 S.E. 51st Street Issaquah, WA 98029, USA" starting from May 11, 2020. The old address of the legal manufacturer can be still indicated in our documents as it was valid before May 11, 2020.

8.3. Additional labeling information for the components;

Name of the component	Label artwork																		
Physio module	Common Physio Module  (01) 0400000005190 GTM (240) 10000000 (422) 276 (21) 1000 (20) 00 MODEL ORIGIN DE SERIAL REV PKK																		
9C3 needle guide	Label Artwork: REF 10136145 LOT XXXXXXX  QTY 1 Biopsy Starter Kit Reusable non-sterile multi-angle bracket with Ultra Pro II™ needle guides (5). (14 x 91.5cm) CIV Files™ covers (5) and cleaning brushes (5) for use with Siemens BC3 HD and 9C3 transducers R ONLY <table><tr><td>(04) Biosoftware</td><td>(11) Kit de soins per biopsie</td><td>(20-CH) 20 個 6.35 x 91.5 ㎝ 針</td></tr><tr><td>(06) Starter biopsie</td><td>(12) 1. 42.2 9.1 5.1 1.1</td><td>(24) 24 個 6.35 ㎝ 針</td></tr><tr><td>(11) Biosoftware</td><td>(40) 40 個 6.35 ㎝ 針</td><td>(54) Biosoftware</td></tr><tr><td>(18) Kit de base pour biopsie</td><td>(60) 60 個 6.35 ㎝ 針</td><td>(78) Biosoftware</td></tr><tr><td>(02) Biosoftware</td><td>(21) 21 個 6.35 ㎝ 針</td><td></td></tr><tr><td>(01) Kit de base pour biopsie</td><td>(00) 00 個 6.35 ㎝ 針</td><td></td></tr></table>  Manufactured for Siemens Medical Systems USA, Inc. 885 East Middlefield Road Mountain View, CA 94031, USA CE 0120 ISO 9001 ISO 13485 ISO 14001 ISO 18001 ISO 26001 ISO 27001 ISO 28001 ISO 29001 ISO 30001 ISO 31001 ISO 32001 ISO 33001 ISO 34001 ISO 35001 ISO 36001 ISO 37001 ISO 38001 ISO 39001 ISO 40001 ISO 41001 ISO 42001 ISO 43001 ISO 44001 ISO 45001 ISO 46001 ISO 47001 ISO 48001 ISO 49001 ISO 50001 ISO 51001 ISO 52001 ISO 53001 ISO 54001 ISO 55001 ISO 56001 ISO 57001 ISO 58001 ISO 59001 ISO 60001 ISO 61001 ISO 62001 ISO 63001 ISO 64001 ISO 65001 ISO 66001 ISO 67001 ISO 68001 ISO 69001 ISO 70001 ISO 71001 ISO 72001 ISO 73001 ISO 74001 ISO 75001 ISO 76001 ISO 77001 ISO 78001 ISO 79001 ISO 80001 ISO 81001 ISO 82001 ISO 83001 ISO 84001 ISO 85001 ISO 86001 ISO 87001 ISO 88001 ISO 89001 ISO 90001 ISO 91001 ISO 92001 ISO 93001 ISO 94001 ISO 95001 ISO 96001 ISO 97001 ISO 98001 ISO 99001 ISO 10001 ISO 10101 ISO 10201 ISO 10301 ISO 10401 ISO 10501 ISO 10601 ISO 10701 ISO 10801 ISO 10901 ISO 11001 ISO 11101 ISO 11201 ISO 11301 ISO 11401 ISO 11501 ISO 11601 ISO 11701 ISO 11801 ISO 11901 ISO 12001 ISO 12101 ISO 12201 ISO 12301 ISO 12401 ISO 12501 ISO 12601 ISO 12701 ISO 12801 ISO 12901 ISO 13001 ISO 13101 ISO 13201 ISO 13301 ISO 13401 ISO 13501 ISO 13601 ISO 13701 ISO 13801 ISO 13901 ISO 14001 ISO 14101 ISO 14201 ISO 14301 ISO 14401 ISO 14501 ISO 14601 ISO 14701 ISO 14801 ISO 14901 ISO 15001 ISO 15101 ISO 15201 ISO 15301 ISO 15401 ISO 15501 ISO 15601 ISO 15701 ISO 15801 ISO 15901 ISO 16001 ISO 16101 ISO 16201 ISO 16301 ISO 16401 ISO 16501 ISO 16601 ISO 16701 ISO 16801 ISO 16901 ISO 17001 ISO 17101 ISO 17201 ISO 17301 ISO 17401 ISO 17501 ISO 17601 ISO 17701 ISO 17801 ISO 17901 ISO 18001 ISO 18101 ISO 18201 ISO 18301 ISO 18401 ISO 18501 ISO 18601 ISO 18701 ISO 18801 ISO 18901 ISO 19001 ISO 19101 ISO 19201 ISO 19301 ISO 19401 ISO 19501 ISO 19601 ISO 19701 ISO 19801 ISO 19901 ISO 20001 ISO 20101 ISO 20201 ISO 20301 ISO 20401 ISO 20501 ISO 20601 ISO 20701 ISO 20801 ISO 20901 ISO 21001 ISO 21101 ISO 21201 ISO 21301 ISO 21401 ISO 21501 ISO 21601 ISO 21701 ISO 21801 ISO 21901 ISO 22001 ISO 22101 ISO 22201 ISO 22301 ISO 22401 ISO 22501 ISO 22601 ISO 22701 ISO 22801 ISO 22901 ISO 23001 ISO 23101 ISO 23201 ISO 23301 ISO 23401 ISO 23501 ISO 23601 ISO 23701 ISO 23801 ISO 23901 ISO 24001 ISO 24101 ISO 24201 ISO 24301 ISO 24401 ISO 24501 ISO 24601 ISO 24701 ISO 24801 ISO 24901 ISO 25001 ISO 25101 ISO 25201 ISO 25301 ISO 25401 ISO 25501 ISO 25601 ISO 25701 ISO 25801 ISO 25901 ISO 26001 ISO 26101 ISO 26201 ISO 26301 ISO 26401 ISO 26501 ISO 26601 ISO 26701 ISO 26801 ISO 26901 ISO 27001 ISO 27101 ISO 27201 ISO 27301 ISO 27401 ISO 27501 ISO 27601 ISO 27701 ISO 27801 ISO 27901 ISO 28001 ISO 28101 ISO 28201 ISO 28301 ISO 28401 ISO 28501 ISO 28601 ISO 28701 ISO 28801 ISO 28901 ISO 29001 ISO 29101 ISO 29201 ISO 29301 ISO 29401 ISO 29501 ISO 29601 ISO 29701 ISO 29801 ISO 29901 ISO 30001 ISO 30101 ISO 30201 ISO 30301 ISO 30401 ISO 30501 ISO 30601 ISO 30701 ISO 30801 ISO 30901 ISO 31001 ISO 31101 ISO 31201 ISO 31301 ISO 31401 ISO 31501 ISO 31601 ISO 31701 ISO 31801 ISO 31901 ISO 32001 ISO 32101 ISO 32201 ISO 32301 ISO 32401 ISO 32501 ISO 32601 ISO 32701 ISO 32801 ISO 32901 ISO 33001 ISO 33101 ISO 33201 ISO 33301 ISO 33401 ISO 33501 ISO 33601 ISO 33701 ISO 33801 ISO 33901 ISO 34001 ISO 34101 ISO 34201 ISO 34301 ISO 34401 ISO 34501 ISO 34601 ISO 34701 ISO 34801 ISO 34901 ISO 35001 ISO 35101 ISO 35201 ISO 35301 ISO 35401 ISO 35501 ISO 35601 ISO 35701 ISO 35801 ISO 35901 ISO 36001 ISO 36101 ISO 36201 ISO 36301 ISO 36401 ISO 36501 ISO 36601 ISO 36701 ISO 36801 ISO 36901 ISO 37001 ISO 37101 ISO 37201 ISO 37301 ISO 37401 ISO 37501 ISO 37601 ISO 37701 ISO 37801 ISO 37901 ISO 38001 ISO 38101 ISO 38201 ISO 38301 ISO 38401 ISO 38501 ISO 38601 ISO 38701 ISO 38801 ISO 38901 ISO 39001 ISO 39101 ISO 39201 ISO 39301 ISO 39401 ISO 39501 ISO 39601 ISO 39701 ISO 39801 ISO 39901 ISO 40001 ISO 40101 ISO 40201 ISO 40301 ISO 40401 ISO 40501 ISO 40601 ISO 40701 ISO 40801 ISO 40901 ISO 41001 ISO 41101 ISO 41201 ISO 41301 ISO 41401 ISO 41501 ISO 41601 ISO 41701 ISO 41801 ISO 41901 ISO 42001 ISO 42101 ISO 42201 ISO 42301 ISO 42401 ISO 42501 ISO 42601 ISO 42701 ISO 42801 ISO 42901 ISO 43001 ISO 43101 ISO 43201 ISO 43301 ISO 43401 ISO 43501 ISO 43601 ISO 43701 ISO 43801 ISO 43901 ISO 44001 ISO 44101 ISO 44201 ISO 44301 ISO 44401 ISO 44501 ISO 44601 ISO 44701 ISO 44801 ISO 44901 ISO 45001 ISO 45101 ISO 45201 ISO 45301 ISO 45401 ISO 45501 ISO 45601 ISO 45701 ISO 45801 ISO 45901 ISO 46001 ISO 46101 ISO 46201 ISO 46301 ISO 46401 ISO 46501 ISO 46601 ISO 46701 ISO 46801 ISO 46901 ISO 47001 ISO 47101 ISO 47201 ISO 47301 ISO 47401 ISO 47501 ISO 47601 ISO 47701 ISO 47801 ISO 47901 ISO 48001 ISO 48101 ISO 48201 ISO 48301 ISO 48401 ISO 48501 ISO 48601 ISO 48701 ISO 48801 ISO 48901 ISO 49001 ISO 49101 ISO 49201 ISO 49301 ISO 49401 ISO 49501 ISO 49601 ISO 49701 ISO 49801 ISO 49901 ISO 50001 ISO 50101 ISO 50201 ISO 50301 ISO 50401 ISO 50501 ISO 50601 ISO 50701 ISO 50801 ISO 50901 ISO 51001 ISO 51101 ISO 51201 ISO 51301 ISO 51401 ISO 51501 ISO 51601 ISO 51701 ISO 51801 ISO 51901 ISO 52001 ISO 52101 ISO 52201 ISO 52301 ISO 52401 ISO 52501 ISO 52601 ISO 52701 ISO 52801 ISO 52901 ISO 53001 ISO 53101 ISO 53201 ISO 53301 ISO 53401 ISO 53501 ISO 53601 ISO 53701 ISO 53801 ISO 53901 ISO 54001 ISO 54101 ISO 54201 ISO 54301 ISO 54401 ISO 54501 ISO 54601 ISO 54701 ISO 54801 ISO 54901 ISO 55001 ISO 55101 ISO 55201 ISO 55301 ISO 55401 ISO 55501 ISO 55601 ISO 55701 ISO 55801 ISO 55901 ISO 56001 ISO 56101 ISO 56201 ISO 56301 ISO 56401 ISO 56501 ISO 56601 ISO 56701 ISO 56801 ISO 56901 ISO 57001 ISO 57101 ISO 57201 ISO 57301 ISO 57401 ISO 57501 ISO 57601 ISO 57701 ISO 57801 ISO 57901 ISO 58001 ISO 58101 ISO 58201 ISO 58301 ISO 58401 ISO 58501 ISO 58601 ISO 58701 ISO 58801 ISO 58901 ISO 59001 ISO 59101 ISO 59201 ISO 59301 ISO 59401 ISO 59501 ISO 59601 ISO 59701 ISO 59801 ISO 59901 ISO 60001 ISO 60101 ISO 60201 ISO 60301 ISO 60401 ISO 60501 ISO 60601 ISO 60701 ISO 60801 ISO 60901 ISO 61001 ISO 61101 ISO 61201 ISO 61301 ISO 61401 ISO 61501 ISO 61601 ISO 61701 ISO 61801 ISO 61901 ISO 62001 ISO 62101 ISO 62201 ISO 62301 ISO 62401 ISO 62501 ISO 62601 ISO 62701 ISO 62801 ISO 62901 ISO 63001 ISO 63101 ISO 63201 ISO 63301 ISO 63401 ISO 63501 ISO 63601 ISO 63701 ISO 63801 ISO 63901 ISO 64001 ISO 64101 ISO 64201 ISO 64301 ISO 64401 ISO 64501 ISO 64601 ISO 64701 ISO 64801 ISO 64901 ISO 65001 ISO 65101 ISO 65201 ISO 65301 ISO 65401 ISO 65501 ISO 65601 ISO 65701 ISO 65801 ISO 65901 ISO 66001 ISO 66101 ISO 66201 ISO 66301 ISO 66401 ISO 66501 ISO 66601 ISO 66701 ISO 66801 ISO 66901 ISO 67001 ISO 67101 ISO 67201 ISO 67301 ISO 67401 ISO 67501 ISO 67601 ISO 67701 ISO 67801 ISO 67901 ISO 68001 ISO 68101 ISO 68201 ISO 68301 ISO 68401 ISO 68501 ISO 68601 ISO 68701 ISO 68801 ISO 68901 ISO 69001 ISO 69101 ISO 69201 ISO 69301 ISO 69401 ISO 69501 ISO 69601 ISO 69701 ISO 69801 ISO 69901 ISO 70001 ISO 70101 ISO 70201 ISO 70301 ISO 70401 ISO 70501 ISO 70601 ISO 70701 ISO 70801 ISO 70901 ISO 71001 ISO 71101 ISO 71201 ISO 71301 ISO 71401 ISO 71501 ISO 71601 ISO 71701 ISO 71801 ISO 71901 ISO 72001 ISO 72101 ISO 72201 ISO 72301 ISO 72401 ISO 72501 ISO 72601 ISO 72701 ISO 72801 ISO 72901 ISO 73001 ISO 73101 ISO 73201 ISO 73301 ISO 73401 ISO 73501 ISO 73601 ISO 73701 ISO 73801 ISO 73901 ISO 74001 ISO 74101 ISO 74201 ISO 74301 ISO 74401 ISO 74501 ISO 74601 ISO 74701 ISO 74801 ISO 74901 ISO 75001 ISO 75101 ISO 75201 ISO 75301 ISO 75401 ISO 75501 ISO 75601 ISO 75701 ISO 75801 ISO 75901 ISO 76001 ISO 76101 ISO 76201 ISO 76301 ISO 76401 ISO 76501 ISO 76601 ISO 76701 ISO 76801 ISO 76901 ISO 77001 ISO 77101 ISO 77201 ISO 77301 ISO 77401 ISO 77501 ISO 77601 ISO 77701 ISO 77801 ISO 77901 ISO 78001 ISO 78101 ISO 78201 ISO 78301 ISO 78401 ISO 78501 ISO 78601 ISO 78701 ISO 78801 ISO 78901 ISO 79001 ISO 79101 ISO 79201 ISO 79301 ISO 79401 ISO 79501 ISO 79601 ISO 79701 ISO 79801 ISO 79901 ISO 80001 ISO 80101 ISO 80201 ISO 80301 ISO 80401 ISO 80501 ISO 80601 ISO 80701 ISO 80801 ISO 80901 ISO 81001 ISO 81101 ISO 81201 ISO 81301 ISO 81401 ISO 81501 ISO 81601 ISO 81701 ISO 81801 ISO 81901 ISO 82001 ISO 82101 ISO 82201 ISO 82301 ISO 82401 ISO 82501 ISO 82601 ISO 82701 ISO 82801 ISO 82901 ISO 83001 ISO 83101 ISO 83201 ISO 83301 ISO 83401 ISO 83501 ISO 83601 ISO 83701 ISO 83801 ISO 83901 ISO 84001 ISO 84101 ISO 84201 ISO 84301 ISO 84401 ISO 84501 ISO 84601 ISO 84701 ISO 84801 ISO 84901 ISO 85001 ISO 85101 ISO 85201 ISO 85301 ISO 85401 ISO 85501 ISO 85601 ISO 85701 ISO 85801 ISO 85901 ISO 86001 ISO 86101 ISO 86201 ISO 86301 ISO 86401 ISO 86501 ISO 86601 ISO 86701 ISO 86801 ISO 86901 ISO 87001 ISO 87101 ISO 87201 ISO 87301 ISO 87401 ISO 87501 ISO 87601 ISO 87701 ISO 87801 ISO 87901 ISO 88001 ISO 88101 ISO 88201 ISO 88301 ISO 88401 ISO 88501 ISO 88601 ISO 88701 ISO 88801 ISO 88901 ISO 89001 ISO 89101 ISO 89201 ISO 89301 ISO 89401 ISO 89501 ISO 89601 ISO 89701 ISO 89801 ISO 89901 ISO 90001 ISO 90101 ISO 90201 ISO 90301 ISO 90401 ISO 90501 ISO 90601 ISO 90701 ISO 90801 ISO 90901 ISO 91001 ISO 91101 ISO 91201 ISO 91301 ISO 91401 ISO 91501 ISO 91601 ISO 91701 ISO 91801 ISO 91901 ISO 92001 ISO 92101 ISO 92201 ISO 92301 ISO 92401 ISO 92501 ISO 92601 ISO 92701 ISO 92801 ISO 92901 ISO 93001 ISO 93101 ISO 93201 ISO 93301 ISO 93401 ISO 93501 ISO 93601 ISO 93701 ISO 93801 ISO 93901 ISO 94001 ISO 94101 ISO 94201 ISO 94301 ISO 94401 ISO 94501 ISO 94601 ISO 94701 ISO 94801 ISO 94901 ISO 95001 ISO 95101 ISO 95201 ISO 95301 ISO 95401 ISO 95501 ISO 95601 ISO 95701 ISO 95801 ISO 95901 ISO 96001 ISO 96101 ISO 96201 ISO 96301 ISO 96401 ISO 96501 ISO 96601 ISO 96701 ISO 96801 ISO 96901 ISO 97001 ISO 97101 ISO 97201 ISO 97301 ISO 9	(04) Biosoftware	(11) Kit de soins per biopsie	(20-CH) 20 個 6.35 x 91.5 ㎝ 針	(06) Starter biopsie	(12) 1. 42.2 9.1 5.1 1.1	(24) 24 個 6.35 ㎝ 針	(11) Biosoftware	(40) 40 個 6.35 ㎝ 針	(54) Biosoftware	(18) Kit de base pour biopsie	(60) 60 個 6.35 ㎝ 針	(78) Biosoftware	(02) Biosoftware	(21) 21 個 6.35 ㎝ 針		(01) Kit de base pour biopsie	(00) 00 個 6.35 ㎝ 針	
(04) Biosoftware	(11) Kit de soins per biopsie	(20-CH) 20 個 6.35 x 91.5 ㎝ 針																	
(06) Starter biopsie	(12) 1. 42.2 9.1 5.1 1.1	(24) 24 個 6.35 ㎝ 針																	
(11) Biosoftware	(40) 40 個 6.35 ㎝ 針	(54) Biosoftware																	
(18) Kit de base pour biopsie	(60) 60 個 6.35 ㎝ 針	(78) Biosoftware																	
(02) Biosoftware	(21) 21 個 6.35 ㎝ 針																		
(01) Kit de base pour biopsie	(00) 00 個 6.35 ㎝ 針																		

**SC1 needle
guide**

Verza™

GUIDANCE SYSTEM

Reusable non-sterile tracking bracket with Verza™ guides (5),
(14 x 91.5cm) CIV-Flax™ covers (5), and cleaning brushes (5)
for use with Siemens Healthineers (Model No. 11291794) SC1 transducers

REF 11287651

LOT XXXXXXX

QTY 1

(8A) Steuerungssystem	(8E) Führungssystem	(8G) 가이드 시스템	(2E-CN) 导引系统
(8B) Cuidado do sistema	(8I) L'entretien du système	(8H) Führungssystem	(8F) Sistema de guía
(8C) Opleidingssysteem	(8J) L'entretien du système	(8I) Führungssystem	(8G) Styringsystem
(8D) Opleidingssysteem	(8K) L'entretien du système	(8J) Führungssystem	(8H) Styringsystem

(8A) Führungssystem	(8E) Führungssystem	(8G) 가이드 시스템	(2E-CN) 导引系统
(8B) Cuidado do sistema	(8I) L'entretien du système	(8H) Führungssystem	(8F) Sistema de guía
(8C) Opleidingssysteem	(8J) L'entretien du système	(8I) Führungssystem	(8G) Styringsystem
(8D) Opleidingssysteem	(8K) L'entretien du système	(8J) Führungssystem	(8H) Styringsystem

SIEMENS HEALTHINEERS SC1

Manufactured for Siemens Medical Solutions USA, Inc.
800 East Main Street, Mountain View, CA 94031 USA

(91) 08841438112277
(11) 0000000
(17) 0000000

Rx ONLY

CE 0120

102 First Street South
Kalama, IA 52247 USA
319.248.8757
WWW.CIVCO.COM

Date of Manufacture
XXXX.XX.XX

Use By
XXXX.XX.XX

11287651-8 REV. C

**DAX needle
guide**

Verza™

GUIDANCE SYSTEM

Reusable non-sterile tracking bracket with Verza™ guides (5),
(14 x 91.5cm) CIV-Flax™ covers (5), and cleaning brushes (5)
for use with Siemens Healthineers DAX transducers

REF 11287652

LOT XXXXXXX

QTY 1

(8A) Steuerungssystem	(8E) Führungssystem	(8G) 가이드 시스템	(2E-CN) 导引系统
(8B) Cuidado do sistema	(8I) L'entretien du système	(8H) Führungssystem	(8F) Sistema de guía
(8C) Opleidingssysteem	(8J) L'entretien du système	(8I) Führungssystem	(8G) Styringsystem
(8D) Opleidingssysteem	(8K) L'entretien du système	(8J) Führungssystem	(8H) Styringsystem

(8A) Führungssystem	(8E) Führungssystem	(8G) 가이드 시스템	(2E-CN) 导引系统
(8B) Cuidado do sistema	(8I) L'entretien du système	(8H) Führungssystem	(8F) Sistema de guía
(8C) Opleidingssysteem	(8J) L'entretien du système	(8I) Führungssystem	(8G) Styringsystem
(8D) Opleidingssysteem	(8K) L'entretien du système	(8J) Führungssystem	(8H) Styringsystem

SIEMENS HEALTHINEERS DAX

Manufactured for Siemens Medical Solutions USA, Inc.
800 East Main Street, Mountain View, CA 94031 USA

(91) 08841438112277
(11) 0000000
(17) 0000000

Rx ONLY

CE 0120

102 First Street South
Kalama, IA 52247 USA
319.248.8757
WWW.CIVCO.COM

Date of Manufacture
XXXX.XX.XX

Use By
XXXX.XX.XX

11287652-8 REV. C

18L6 needle
guide

Label Artwork:

REF 10044200 **LOT XXXXXXX**

QTY 1

Biopsy Starter Kit

Reusable non-sterile multi-angle bracket with Ultra Pro II™ needle guides (5)
(14 x 9) 5cm CIV-Flar™ covers (5) and cleaning brushes (5) for use with
Siemens 18L6 HD transducers

R ONLY

(04) Biopsy-startpakete	(17) Kit di avvio per biopsia	(28-CN) 活組織切片検査套件
(05) Starter biopsie	(18) 开始スターキット	(29) Equipo inicial para biopsia
(06) Biopsy-startpakete.s	(20) 开始 스타터 키트	(30) Biopsy Starter Kit
(07) Kit de base pour biopsies	(21) Starter for biopsy	(31) Biopsy Starter Kit
(08) Biopsy-Startset	(22) Kit de biopsia	
(09) Kit Biopsie Starter	(23) Transducerul necesare pentru biopsie	

Manufactured for
Siemens Medical Solutions USA, Inc.
605 East Middlefield Road
Mountain View, CA 94043, USA

10100947 4301 00238
10100000000
11100000000
11100000000

CE 0120

Temp. range
-10°C to +40°C

Keep away from sunlight

Consult instructions for use

CIVCO
102 First Street South
Kalamazoo, MI 49001 USA
319 248 8757
WWW.CIVCO.COM

10044200 REV E

14L5 needle
guide

Label Artwork:

Verza
GUIDANCE SYSTEM

REF 11287650
LOT XXXXXXX
QTY 1

Reusable non-sterile bracket with Verza™ guides (5) and
(14 x 9) 5cm CIV-Flar™ covers (5) for use with
Siemens Healthineers 14L5 (Model No. 11254C90) transducers

(04) VerzaSystem	(05) VerzaSystem	(06) VerzaSystem	(07) VerzaSystem
(08) VerzaSystem	(09) VerzaSystem	(10) VerzaSystem	(11) VerzaSystem
(12) VerzaSystem	(13) VerzaSystem	(14) VerzaSystem	(15) VerzaSystem
(16) VerzaSystem	(17) VerzaSystem	(18) VerzaSystem	(19) VerzaSystem
(20) VerzaSystem	(21) VerzaSystem	(22) VerzaSystem	(23) VerzaSystem
(24) VerzaSystem	(25) VerzaSystem	(26) VerzaSystem	(27) VerzaSystem

CIVCO

SIEMENS HEALTHINEERS 14L5

Manufactured for Siemens Medical Solutions USA, Inc.
605 East Middlefield Road, Mountain View, CA 94043, USA

10100947 4301 00238
10100000000
11100000000
11100000000

CE 0120

Temp. range
-10°C to +40°C

Keep away from sunlight

Consult instructions for use

CIVCO
102 First Street South
Kalamazoo, MI 49001 USA
319 248 8757
WWW.CIVCO.COM

11287650 REV C

10L4 needle
guide

Label Artwork:

[illegible][illegible]

Label Artwork:

LOT XXXXXXX

CIVCO

Needle Guide

Assembly non sterile needle guide with cleaning brushes / 55
for use with Samara OFC, MAC3 and EC9 for Introduttore

2, ONLY

[BA] Nadelguide	[ET] Guida di ago	[ZH-CN] 針導引器
[NL] Naaldleider	[JA] ナイロガイド	[ES] Guía de aguja
[FI] Neulanjohti	[KO] 바늘 가이드	[SV] Nålguide
[FR] Ode de-aiguille	[NO] Nålelydning	[IT] Igne-Aiguiz
[DE] Nadelführung	[PT-PT] Guia de agulhas	
[EL] Οδηγός Βελόνας	[RU] Трассовод иглы	

Manufactured for
SARAFIM GROUP, Sarafim Inc. via
CIVCO Medical Products (Pvt)
Houston, Texas, USA 77058, USA

CE
0126

RoHS
XAUCE-001X

CIVCO
100 River Street South
Houston, TX 77058 USA
281 348 0707
WWW.CIVCO.COM

040000014000000A
P1100000000000
11000000000000

NOT FOR SALE IN THE UNITED STATES OF AMERICA

©2006 CIVCO INC. ALL RIGHTS RESERVED

4V1 needle
guide

Label Artwork:

REF 10441187

QTY 1

LOT XXXXXXX

Tracking Bracket Starter Kit
Reusable non-sterile multi-angle tracking bracket with Ultra-Pro II™ needle guides (5), (14 x 97.5cm) CN Flex™ covers (5), cleaning brushes (5) for use with Siemens 4V1 transducers

R ONLY

(04) Springklemm-Startsat	(11) Kit di avvio per immobilizzazione paziente	(24) (K) 患者固定具
(05) Startkit voor opsporingsteugel	(34) トラッキングブラケットスターキット	(35) Equipo básico de soporte de cartita
(11) Anslutningspaket	(40) 患者固定具セット	(51) Startkits for opsporingsteugel
(18) Kit de démarrage de la console de surveillance	(40) Startkit for opsporingsteugel	(78) Starter Cesteği Başlangıç Kit
(06) Startset Tracking-uitrusting	(91-97) Kit de suporte de rastreo	
(41) Kit d'attachement, immobilisateur, Déplacement	(98) Handreuk voor het openen van de tracking-uitrusting	

Manufactured for
Siemens Medical Solutions USA, Inc.
999 Egan Highway Road
Mountain View, CA 94035 USA

CE 0120

Keep away from sunlight

Consult instructions for use

Date of Manufacture
XXXX XX XX

Expiry Date
XXXX XX XX

CIVCO
102 First Street South
Kalamazoo, MI 49001 USA
313 248 8737
WWW.CIVCO.COM

Footswitch

KINESIS
COMPUTER ERGONOMICS
Phone: 800 454 6374 www.kinesis.com

Savant™ Elite2™
Waterproof Triple Foot Pedal

PN: FPW30A-IPX8-S2
SN: 30024UPW3-S2

RoHS IPX8

FC Tested To Comply With FCC Standards FOR HOME OR OFFICE USE

Assembled in USA by Kinesis® Corporation

CE

5V DC 50mA

Intertek 4007055

9EC4 needle
guide,
disposable

Label Artwork:

REF 10659482

QTY 24

LOT XXXXXXX

Needle Guide

Sterile endocavity needle guide with (2.6 x 20cm) and (3.5 x 20cm) latex covers
for use with Siemens 9EC4, 9MC3, EC9 4w, and MC9 4 transducers

[de] Nadelguide
[nl] Naaldvoert
[fi] Naidurjoht
[es] Guías de aguja
[pt] Guia do agulhas
[fr] Cône de guidage
[it] Guida dell'ago
[ru] Ланцет-руководящий
[uk] Ланцет-руководящий
[pl] Lancyj przewodnik
[tr] İğne Kılavuzu
[zh-CN] 针导引器
[ja] 針導引器
[ko] 바늘 안내기
[he] מדריך לכף
[ar] دليل الإبرة
[id] Petunjuk jarum
[vi] Hướng dẫn kim
[th] เข็มนำทาง
[he] מדריך לכף
[ar] دليل الإبرة
[id] Petunjuk jarum
[vi] Hướng dẫn kim
[th] เข็มนำทาง

[de] Nadelguide
[nl] Naaldvoert
[fi] Naidurjoht
[es] Guías de aguja
[pt] Guia do agulhas
[fr] Cône de guidage
[it] Guida dell'ago
[ru] Ланцет-руководящий
[uk] Ланцет-руководящий
[pl] Lancyj przewodnik
[tr] İğne Kılavuzu
[zh-CN] 针导引器
[ja] 針導引器
[ko] 바늘 안내기
[he] מדריך לכף
[ar] دليل الإبرة
[id] Petunjuk jarum
[vi] Hướng dẫn kim
[th] เข็มนำทาง

WARNING:
Components of this product contain NATURAL RUBBER (LATEX) which may cause allergic reactions in some individuals.

Manufactured for:
Siemens Medical Solutions USA, Inc.
880 East Mainfield Road
Mountain View, CA 94043, USA

CE 0120

Do not use if damaged or deformed

CE 0120

Do not use if damaged or deformed

CE 0120

Do not use if damaged or deformed

DATE OF MANUFACTURE: XXXXX.XX.XX
USE BY: XXXXX.XX.XX

CIVCO
102 First Street South
Kalamazoo, MI 49001 USA
319.248.6757
WWW.CIVCO.COM

10659482

eSie Fusion
accessories

(01)0000000000000000 GTIN

(20)00 Rev

(21)000000 Serial

(240)00000000 Model

(422)840 Origin US

IVK

SIEMENS

(1P) Model No. 10441625

(S) Serial No. 3283

(2P) Rev. 02

IVK

ECG cable
(EU type)

Speciality for getting
and transferring signals

Device name: ECG CABLE WITH LEADS

Model No.: 11148104

REF: EA023C5I

LOT:

BN: See Cable

Quantity: 1 Piece

Compatible: SIEMENS P/N 11148104 REV 05

 Shenzhen Med-link Electronics Tech Co., Ltd.
4th Floor, Building A, Ying Tailing Industrial Park,
Dalong South Road, Longhua District, Shenzhen,
Guangdong, 518109, CHINA
<http://www.med-linket.com>
E-mail: sales@med-linket.com
Made in China

 Shanghai International Hoisting Corp. GmbH (Europe)
Erfenstrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175 Fax: +49-40-255728


9. Additional Storage information

Needle guides for 9C3, 5C1, DAX, 18L6, 14L5, 10L4, 4V1 transducers

	During Storage or Transportation
Temperature:	+11°C to +29°C
Atmospheric pressure:	No specific requirements
Relative humidity:	No specific requirements

Disposable needle guides for 9EC4 transducers

	During Storage or Transportation
Temperature:	+11°C to +29°C
Atmospheric pressure:	No specific requirements
Relative humidity:	No specific requirements

Reusable needle guides for 9EC4 transducers

	During Storage or Transportation
Temperature:	+11°C to +29°C
Atmospheric pressure:	No specific requirements

Relative humidity:

No specific requirements

10. Service Life

The expected service life of the system is 10 years

11. Additional maintenance information.

Customer's support center of Siemens Healthcare LLC:

Tel: 8 (800) 550 5022 (24 hours call center)

Email: service.ru@siemens-healthineers.com

12. Warranty

The factory warranty is provided for the period of 12 months after commissioning of the equipment but not more than 15 months from the data of shipment by the manufacturer.


13. Complaints

In the Russian Federation, for questions regarding complaints, repair or maintenance of the ACUSON Sequoia Diagnostic Ultrasound System with Accessories manufactured by Siemens Medical Solutions USA, Inc., USA, contact Siemens Healthcare LLC, legal and post address: Dubininskaja street 96, 115093, Moscow, Russia, Phone: (495) 737-12-52; Fax: (495) 737-13-20.

Glossary

Name of the component in application form for Russian registration	Name of the component in IFU
Needle guide bracket for 10L4 transducer	Biopsy and Fusion needle guide set, 10L4
Needle guide bracket for 5C1 transducer	Biopsy and Fusion needle guide set, 5C1
Needle guide bracket for DAX transducer	Biopsy and Fusion needle guide set, DAX
Needle guide bracket for 14L5 transducer	Biopsy needle guide set, 14L5
Needle guide bracket for 4V1 transducer	Fusion needle guide set, 4V1
Needle guide bracket for 18L6 transducer	Biopsy needle guide set, 18L6
Needle guide bracket for 9C3 transducer	Biopsy needle guide set, 9C3
Reusable Needle guide bracket for 9EC4 transducer	Reusable endocavity needle guide, 9EC4
Disposable Needle guide bracket for 9EC4 transducer	Disposable endocavity needle guide, 9EC4
eTrax sensors	CIVCO General Purpose 8 mm sensors
Module for cardiac examinations with function of cardiologic measures and reports	Cardiac imaging option

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Verza™

GUIDANCE SYSTEM

Reusable non-sterile tracking bracket with Verza™ guides (5),
 (14 x 9) 5cm CIV-Flux™ covers (5), and cleaning brushes (5)
 for use with Siemens Healthineers DAX transducers.

(0A) Stylsystem	(0E) Führungssystem	(0B) 가이드 시스템	(0D) Reinigungssystem	(ZB-CB) 가이드 시스템	(SS) Sistema de guía	(TS) Kitverba System
(0L) Oeko dringssystem	(0I) Leuchtgeleitsystem	(NB) Führungssystem	(ST-ST) Sistema de limpieza	(SB) Styrmagnetsystem		
(FJ) Ökoreinigungssystem	(IT) Sistema guida	(RW) Cactores para desengrasar actua				
(FR) Système de guidage	(JA) ガイドシステム					

SIEMENS HEALTHINEERS DAX

Manufactured for Siemens Medical Solutions USA, Inc.
 One East Wacker Drive, Deerfield, IL 60015, USA

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QTY 1

Biopsy Starter Kit

Reusable non-sterile multi-angle bracket with Ultra Pro II™ needle guides (5)
(14 x 21.5 cm), CIV Flex™ covers (5), and cleaning brush (5), for use with
Siemens 18L6 HD transducers

R ONLY

[DA] Biopsy starter kit	[ET] Kit de inicio per biopsia	[ZH-CN] 活检开始套件
[NL] Starter biopsie	[JA] 生検開始キット	[ES] Equipo inicial para biopsia
[FI] Biopsia-alku, aloittelus	[KO] 생검 스타터 키트	[SV] Biopsi Hyrtärar Kit
[FR] Kit de base pour biopsies	[NO] startersett for biopsi	[TR] Biopsi Başlangıç Kit
[DE] Biopsie-Startset	[PT-PT] Kit de biopsia	
[EL] Kit Εξοικονόμησης Βιοψίας	[RU] Начальный набор для биопсии	

Manufactured for
Siemens Medical Solutions USA, Inc.
805 East Middlefield Road
Mountain View, CA 94043, USA

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GUIDANCE SYSTEM

REF 11287650 **LOT XXXXXXXX** **QTY 1**

Reusable non-sterile bracket with Verza™ guides (5) and
(14 x 21.5 cm), CIV Flex™ covers (5), for use with
Siemens Healthineers 14L5 (Model No. 11254C90) transducers

[DA] Sysystem	[DE] Führungssystem	[KO] 시스템	[ZH-CN] 引导系统
[NL] Leidingssysteem	[FI] Johtosysteemi	[JA] 導引システム	[ES] Sistema de guía
[FR] Outils système	[IT] Sistema guida	[PT-PT] Sistema de orientação	[SV] Styningssystem
[FR] Système de guidage	[JA] 導引システム	[RU] Система для введения иглы	[TR] Kılavuz Sistemi

CIVCO

SIEMENS HEALTHINEERS 14L5

Manufactured for Siemens Medical Solutions USA, Inc.
805 East Middlefield Road, Mountain View, CA 94043, USA

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Verza™

GUIDANCE SYSTEM

Reusable non-sterile tracking bracket with Verza™ guide (1)
(14 x 9" Slim, CIV Hse™ system (2), and cleaning brushes (3)
for use with Siemens Healthineers 10L4 transducers

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QTY 1

SIEMENS HEALTHINEERS 10L4

Manufactured by Siemens Medical Systems USA, Inc.
100 First Street South
Kalamazoo, MI 49001 USA

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Label Artwork:

Needle Guide

Reusable non-sterile needle guide with cleaning brushes (3)
for use with Siemens 9EC4, 9MC2 and EC9 4w transducers

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CE
5V DC 50mA

RECOGNIZED COMPONENT

FCC Tested To Comply
With FCC Standards
FOR HOME OR OFFICE USE

ETL US
C

Assembled in USA by
Kinesis® Corporation

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REF 10659482 LOT XXXXXXXX
QTY 24

Needle Guide

Sterile and/or single needle guide with (2.0 x 25mm) and (2.0 x 25mm) other sizes
For use with Siemens 9EC4, 9MC3, ECR 4w, and MC3 4 transducers

STERILE EO
Sterilization using ethylene oxide
R ONLY

(0A) Needle guide	(11) Guide deflating	(20) (10) (11) (12)
(0B) Needle guide	(14) (11) (12) (13)	(25) Guide deflating
(0C) Needle guide	(15) (11) (12) (13)	(30) Needle
(0D) Needle guide	(16) Needle	(35) Needle
(0E) Needle guide	(17) (11) (12) (13) (14)	(40) Needle
(0F) Needle guide	(18) Needle	(45) Needle

WARNING: Carefully follow the product label and the instructions for use. Do not use the product after the expiration date or if the package is damaged.

Manufactured by
Siemens Medical Solutions USA, Inc.
400 East Main Street
Madison, WI 53703, USA

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
		Speciality for getting and transferring signals	
Device name: ECG CABLE WITH LEADS			
Model No.: 11148104			
REF	EA023C5I		
LOT			
EN	See Cable		
Quantity: 1 Piece			
Compatible: SIEMENS P/N 11148104 REV 05			
 Shenzhen Med-link Electronics Tech Co., Ltd. 4th Floor, Building A, Ying Tallong Industrial Park, Dalong South Road, Longhua District, Shenzhen, Guangdong, 518109, CHINA http://www.med-linket.com E-mail: sales@med-linket.com Made in China			
  Shanghai International Holding Corp. GmbH (Europe) Elffootssee 80, 20537 Hamburg, Germany Tel: +49-40-2513175 Fax: +49-40-255728			
			
			

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AM

Commonwealth of Virginia
The Honorable Kelly Thomasson - Secretary of the Commonwealth

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Aida A. Naanou, certify that the attached document which is a Self-Truthfulness Declaration for Product Registration "Acuson Sequoia" received from the company Siemens Medical Solutions, USA Inc. signed by Hyunjung Lee is an original document.

Aida A Naanou

City/County of Fairfax
Commonwealth of Virginia
Subscribed and sworn to before me
this 14th day of May, 2020
by Aida A Naanou
Perla Kara Al-Shayeb Notary Public
PERLA KARA AL-SHAYEB
Reg # 7338574 Commission expires 1/31/20 22

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- *Advanced Imaging Manual for ACUSON Sequoia system.*

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
- *Instruction for use for needle guides for 5C1 transducer*

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Sincerely yours,

Siemens Medical Solutions USA, Inc.

Hyunjung Lee
Regulatory Affairs Specialist
Siemens Medical Solutions USA, Inc.

May 1, 2020

Date

Verza™


GUIDANCE SYSTEM

REFERENCE GUIDE

Tracking Bracket
for use with
Siemens Healthineers
5C1
transducers

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


IMAGING GUIDANCE SYSTEM


1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

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
ECTING AND STERILIZING BRACKET

Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedures.

	Country of Origin	Manufacturer	Qualified Use	Active Ingredient	Disinfectant Type
X® or equivalent 2% aldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Glutaraldehyde	HLD, S
EX® OPA or equivalent .55% -phthalaldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Ortho-phthalaldehyde	HLD
®	US	ASP® Advanced Sterilization Products	Pre-cleaner	Enzymes, surfactant	CL
sepi® FF (new)	DE	Schülke	Soak	Succinic dialdehyde	HLD

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


LEANING THE GENERAL PURPOSE SENSOR

1. After each use remove general purpose sensor, discarding single-use, disposable components. Remove visible contaminants from general purpose sensor surface.
2. Inspect general purpose sensor for damage such as cracks or breakage.

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OTE: Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedures.

	Country of Origin	Manufacturer	Qualified Use	Active Ingredient	Disinfectant Type
DEX® or equivalent 2% glutaraldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Glutaraldehyde	HLD, S
DEX® OPA or equivalent .55% ortho-phthalaldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Ortho-phthalaldehyde	HLD
zol®	US	ASP® Advanced Sterilization Products	Pre-cleaner	Enzymes, surfactant	CL

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[Signature]

Commonwealth of Virginia
The Honorable Kelly Thomasson - Secretary of the Commonwealth
Document revised April 18, 2016

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This Apostille is not valid for use anywhere within the United States of America or its territories.

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
Notary Public
PERLA KARA AL-SHAYER
Reg # 7338574 Commission expires 1/31/20

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
- *Instruction for use for needle guides for 9C3, 14L5 and 18L6 transducers*

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




BRACKET / NEEDLE GUIDE

INDICATIONS FOR USE

The guidance system is intended for directing instruments such as catheters, electrodes and needles into a targeted anatomical location of a patient relative to the imaging instrument for percutaneous procedures. Prior to use healthcare workers should be trained in ultrasonography.

INTENDED USE

The bracket and disposable guide provide physicians with a tool for performing needle-guided (or catheter) procedures with the use of diagnostic ultrasound transducers.

COVER

INDICATIONS FOR USE / INTENDED USE

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides means for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

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WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- CIVCO's needle guidance systems are designed and validated as an integrated two-part system consisting of a custom reusable bracket and a disposable needle guide. Use of non-CIVCO guides on a CIVCO bracket may result in serious injury to the patient and/or user. To ensure safe and effective use in compliance with system software guidelines, only CIVCO needle guides are recommended.
- Bracket is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned, sterilized or disinfected before each use.
- Do not use if bracket is damaged or does not fit properly.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- For illustration purposes only, transducer and bracket may be shown without a transducer cover. Always place a cover over transducer and bracket to protect patients and users from cross-contamination.

NOTE: The use of this needle guide is not necessarily the same as needle guides for other Siemens transducers. For use of this product, please follow carefully the instructions in this reference guide.

NEEDLE PATH VERIFICATION


- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system user guide.
- For multi-angle brackets, verify bracket angle corresponds with system software before each use.

NOTE: Use appropriate needle length to reach target area.

USING NEEDLE GUIDANCE SYSTEM


1. Attach bracket to transducer by aligning locating features. Ensure bracket is firmly attached. (See illustrations on page 2.)

- 2. For multi-angle brackets only.
 - I. Pull up lock pin.
 - II. Move to desired angle.
 - III. Press lock pin securely into position.


(For A/B references, see table on page 1.)

- 3. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
- 4. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
- 5. Secure cover with enclosed bands.
- 6. Inspect cover to ensure there are no holes or tears.


- 7. Using proper sterile technique, snap unlocked needle guide onto attachment area of bracket.
- 8. Push lock into locked position.
- 9. Select appropriate size needle guide insert and slide into position.

QUICK RELEASE FUNCTION

- 1. Press tab toward bracket to activate needle guide quick-release function.

CLEANING, DISINFECTING AND STERILIZING BRACKET


WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.
- See your system's user guide for reprocessing transducer between use.
- This procedure for reprocessing has been validated for effectiveness and compatibility. Device could be damaged or cross-contaminated due to improper reprocessing.

- Do not gas sterilize or autoclave bracket.
- Inspect bracket for damage such as cracks or breakage. If damage is evident, discontinue use of bracket and contact your ultrasound system support representative.

ASSEMBLING BRACKET FOR CLEANING

NOTE: For multi-angle brackets only.

CLEANING THE BRACKET

1. After each use, remove bracket and needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
2. Soak bracket for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL® Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove bracket from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

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- High-level disinfect bracket using CIDEX OPA® ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize bracket using CIDEX® Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus® (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.

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CIVCO

CORPORATE HEADQUARTERS:
2301 Jones Blvd., Carrollville, IA 52241
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

USA OFFICE:
102 First Street South, Kolono, IA 52247-9589
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

COPYRIGHT © 2018 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ULTRA-PRO IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

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Aida A Naanou

City/County of Fairfax
Commonwealth of Virginia
Subscribed and sworn to before me
this 14th day of May, 2020
by Aida A Naanou
Perla Kara Al-Shayeb Notary Public
PERLA KARA AL-SHAYEB
Reg # 7338574 Commission expires 1/31/20

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
- *Instruction for use for disposable needle guides for 9EC4 transducers*

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USABLE ENDOCAVITY NEEDLE GUIDE

REFERENCE GUIDE

For use with
Siemens
9EC4 / 9MC3 / EC9-4w / MC9-4
transducers

CE
0120

JSE

Needle / instrument guided procedures - transvaginal and transrectal: tissue biopsy, fluid aspiration, and treatment.

A mechanical means for performing needle / instrument guided procedures with the use of the endocavity transducer. Needle guide is attached over endocavity transducer / probe / scanhead and provides a fixed path for the needle or instrument that when coupled by the ultrasound system provides imaging guidelines for visualizing guided instrument placement procedures. CIVCO Endocavity Biopsy / Guides are furnished sterile; single use patient / procedure, disposable. The single use, helps to prevent transfer of microorganisms, body fluids, and particle material to the patient and during reuse of the transducer.

FOR USE / INTENDED USE

Sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows easier scanning and needle guided procedures for body surface, endocavity, and intra-operative use, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides maintenance of a sterile field (sterile covers only). CIVCO Latex Transducer Covers are furnished sterile & single use patient / procedure, disposable.

(United States) law restricts this device to sale by or on the order of a physician.

For use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your transducer's user guide.

Transducer covers containing natural rubber latex may cause allergic reactions. In the USA, refer to FDA Medical Alert MDA91-1.

Use of a needle guide on an endocavity transducer may increase the risk of vaginal or rectal abrasion.

For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.

Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.

Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.

Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.

Do not use if needle guide is damaged or does not fit properly.

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See transducer user guide for proper procedure.

Use care when inserting needle to avoid bending.

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
NEEDLE GUIDANCE

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.


4. Inspect cover to ensure there are no holes or tears.


5. Using proper sterile technique, align guide with locating features of transducer.

6. Press down to secure to transducer.

7. If desired, apply secondary cover by following the steps above.

DISPOSAL

- Dispose of single-use components as infectious waste.
- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's user guide for reprocessing transducer between use.

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
- *Instruction for use for reusable needle guides for 9EC4 transducers*

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REUSABLE ENDOCAVITY NEEDLE GUIDE

REFERENCE GUIDE

For use with
Siemens
9EC4 / 9MC3 / EC9-4w
transducers

CE
0120

NEEDLE GUIDE**INDICATIONS FOR USE / INTENDED USE**

This device provides a mechanical means for performing needle / instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer.

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WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Use of a needle guide on an endocavity transducer may increase the risk of vaginal or rectal abrasion.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Needle guide is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure needle guide is properly cleaned, sterilized or disinfected before each use.
- Do not use if needle guide is damaged or does not fit properly.

NEEDLE PATH VERIFICATION

Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system user guide.

E: Use appropriate needle length to reach target area.

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



With clamp open, align guide with locating features on transducer.

6. Close needle guide around transducer and check for secure fit.

CLEANING, DISINFECTING AND STERILIZING NEEDLE GUIDE

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.
- See your system's user guide for reprocessing transducer between use.
- This procedure for reprocessing has been validated for effectiveness and compatibility. Device could be damaged or cross-contaminated due to improper reprocessing.
- Inspect needle guide for damage such as cracks or breakage. If damage is evident, discontinue use of needle guide and contact your ultrasound system support representative.
- Due to the size and length of the needle guide, an automated method for cleaning and disinfecting cannot be specified.

CLEANING THE NEEDLE GUIDE

1. After each use, remove needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from needle guide surface using a small, soft instrument brush. Keep needle guide from drying out until complete cleaning can be accomplished.
2. Soak needle guide for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL® Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove needle guide from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

DISINFECTING AND STERILIZING NEEDLE GUIDE

- High-level disinfect needle guide using CIDEX OPA® ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehydebased solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize needle guide using CIDEX® Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus® (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Needle guide may be sterilized by using steam (autoclave) at 132-138°C for 15 to 30 minutes.

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Aida A Naanou

City/County of Fairfax
Commonwealth of Virginia
Subscribed and sworn to before me
this 14th day of May, 2020
by Aida A Naanou
[Signature] Notary Public
PERLA KARA AL-SHAYEB
Reg.# 7338574 Commission expires 1/31/2022

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
- *Instruction for use for needle guides for 4V1 transducers*

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




CKET / NEEDLE GUIDE

ICATIONS FOR USE

guidance system is intended for directing instruments such as catheters, electrodes and needles into a targeted tomical location of a patient relative to the imaging instrument for percutaneous procedures. Prior to use healthcare ers should be trained in ultrasonography.

ENDED USE

bracket and disposable guide provide physicians with a tool for performing needle-guided (or catheter) procedures the use of diagnostic ultrasound transducers.

VER

ICATIONS FOR USE / INTENDED USE

active cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative nstic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the ent and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides eans for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished ile & non-sterile; single use patient / procedure, disposable.

UTION

eral (United States) law restricts this device to sale by or on the order of a physician.

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
2. For multi-angle brackets only.
- I. Pull up lock pin.


II. Move to desired angle.

III. Press lock pin securely into position.


(For A/B references, see table on page 1.)

3. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
4. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
5. Secure cover with enclosed bands.
6. Inspect cover to ensure there are no holes or tears.


7. Using proper sterile technique, snap unlocked needle guide onto attachment area of bracket.
8. Push lock into locked position.
9. Select appropriate size needle guide insert and slide into position.

QUICK RELEASE FUNCTION

1. Press tab toward bracket to activate needle guide quick-release function.

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
NOTE: For multi-angle brackets only.

CLEANING THE BRACKET

- After each use, remove bracket and needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
- Soak bracket for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL® Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove bracket from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

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
, Aida A. Naanou, certify that the attached document which is a Self-Truthfulness Declaration for Product Registration " Acuson Sequoia" received from the company Siemens Medical Solutions, USA Inc. signed by Hyunjung Lee is an original document.

Aida A Naanou

City/County of Fairfax
Commonwealth of Virginia
Subscribed and sworn to before me
this 14th day of May, 2020
by Aida A Naanou

Notary Public
PERLA KARA AL-SHAYEB
Reg.# 7338574 Commission expires 1/31/2022

Self – Truthfulness Declaration for Product Registration

Manufacturer: Siemens Medical Solutions USA, Inc.

Product: ACUSON Sequoia

As responsible individuals authorized to represent and act on behalf of Siemens Medical Solutions USA, Inc., we hereby certify that the document:

- *Instruction for use for needle guides for 10L4 transducers*


is genuine, complete and valid.

We shall bear responsibility arising from false, incomplete or invalid documents.

Sincerely yours,
Siemens Medical Solutions USA, Inc.

HyunJung Lee
Regulatory Affairs Specialist
Siemens Medical Solutions USA, Inc.

Date

Verza™




GUIDANCE SYSTEM

REFERENCE GUIDE

Tracking Bracket
for use with
Siemens Healthineers
10L4
transducers

CE
0120

TRACKING BRACKET**INDICATIONS FOR USE / INTENDED USE**

The device is intended to provide physicians with tools for electromagnetic tracking of instruments with respect to image data.

GUIDANCE SYSTEM**INDICATIONS FOR USE / INTENDED USE**

The system provides guidance for precise instrument placement of common interventional devices by positioning the device relative to the ultrasound transducer and the resulting image during a diagnostic or therapeutic procedure. This guidance system is intended for use with pediatric and adult patients.

COVER**INDICATIONS FOR USE / INTENDED USE**

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Verza Guidance System is designed and validated for use with transducer(s) listed on the cover of this instruction.
- Tracking bracket is for use with CIVCO general purpose sensors only.
- General purpose sensor is for use with CIVCO accessories only.
- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your ultrasound system's user guide.
- Prior to use, ensure general purpose sensor is tracking properly.
- To avoid possible position distortion, keep general purpose sensor away from metal.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- CIVCO's guidance systems are designed and validated as an integrated two-part system consisting of a custom reusable bracket and a disposable guide. Use of non-CIVCO guides on a CIVCO bracket or locating feature may result in serious injury to the patient and/or user. To ensure safe and effective use in compliance with system software guidelines, only CIVCO guides are recommended.
- Bracket and general purpose sensor are packaged non-sterile and are reusable. To avoid possible patient contamination, ensure bracket and general purpose sensor are properly cleaned and sterilized before each use.
- Do not use if bracket or general purpose sensor appear damaged or do not fit properly.
- Do not use if locating feature is damaged or guide does not fit properly.
- Ensure gauge insert is in closed position prior to beginning procedure.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- For illustration purposes only, components may be shown without a cover. Always place a cover over transducer, general purpose sensor and bracket to protect patients and users from cross-contamination.

NEEDLE PATH VERIFICATION

- Before using component for first time, perform needle path verification to verify system and Verza guide relationships, as described in ultrasound system user guide.
- Before each use, verify guidance system angles (5) correspond with system software.


NOTE: Use appropriate needle length to reach target area. Verza guide channel measures 3.4 cm (1.34 in).

WALKING TRACKING SYSTEM

- With latch open, align bracket with locating features on transducer.
- Close latch around transducer and click to secure. Ensure bracket is securely attached.

- Place general purpose sensor onto sensor mount.
- Press general purpose sensor toward bracket and click to secure.


WARNING


Ensure general purpose sensor is securely attached. Failure to do so could result in inaccuracy.

WALKING GUIDANCE SYSTEM


- Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
- Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
- Secure cover with enclosed bands.
- Inspect cover to ensure there are no holes or tears.

- Remove Verza from guide storage column of gauge insert tree.

- Using proper sterile technique, place Verza guide onto VerzaLink locating feature.

- Rotate guide locking lever into locked position. Ensure guide is securely attached.

8. FIXED ANGLE APPROACH:


Using angle selector grip, slide angle indicator window to desired position and push angle selector lock toward transducer to secure.

NOTE: Ensure angle selector lock is secure in desired position.

FREE ANGLE APPROACH:

Leave angle selector lock in unlocked position and instrument remains in-plane within constraints of guide angle range.

9. Select appropriate size gauge insert and slide into guide body chute.

BOLT-ACTION QUICK-RELEASE

1. Rotate quick-release tab to open guide insert funnel.

REMOVING GUIDANCE SYSTEM

1. Rotate guide locking lever to unlocked position to remove guide from VerzaLink locating feature of bracket.

2. Remove and dispose of single-use components according to facility infection control policies.
3. Lift bracket latch to remove from transducer. Retain bracket for reprocessing.

REMOVING GENERAL PURPOSE SENSOR FROM BRACKET

1. Pull general purpose sensor away from bracket to remove.

CLEANING, DISINFECTING AND STERILIZING

WARNING

- Users of these products have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single-use components as infectious waste. Clean reusable components after each use.
- See your ultrasound system's user guide for reprocessing transducer between use.
- These procedures for reprocessing have been validated for effectiveness and compatibility. Devices could be damaged or cross-contaminated due to improper reprocessing.
- Do not gas sterilize or autoclave bracket or general purpose sensor.
- If bracket or general purpose sensor damage is evident, discontinue use and contact your CIVCO representative.

CLEANING THE BRACKET

- After each use, remove bracket from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
- Inspect bracket for damage such as cracks or breakage.

DISINFECTING AND STERILIZING BRACKET

- TE: Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedures.

Disinfection	Country of Origin	Manufacturer	Qualified Use	Active Ingredient	Disinfectant Type
DEX® or equivalent 2% glutaraldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Glutaraldehyde	HLD, S
DEX® OPA or equivalent .55% ortho-phthalaldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Ortho-phthalaldehyde	HLD
Enzymatic Cleaner®	US	ASP® Advanced Sterilization Products	Pre-cleaner	Enzymes, surfactant	CL
Sept® FF (new)	DE	Schülke	Soak	Succinic dialdehyde	HLD

DISINFECTANTS AND CLEANING SOLUTIONS TABLE LEGEND

Country of Origin	Disinfectant Type
DE = Germany	CL = Cleaner
US = United States	HLD = High-level disinfectant
	S = Sterilant




CLEANING THE GENERAL PURPOSE SENSOR

- After each use remove general purpose sensor, discarding single-use, disposable components. Remove visible contaminants from general purpose sensor surface.
- Inspect general purpose sensor for damage such as cracks or breakage.

DISINFECTING AND STERILIZING GENERAL PURPOSE SENSOR

- TE: Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedures.

Disinfection	Country of Origin	Manufacturer	Qualified Use	Active Ingredient	Disinfectant Type
DEX® or equivalent 2% glutaraldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Glutaraldehyde	HLD, S
DEX® OPA or equivalent .55% ortho-phthalaldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Ortho-phthalaldehyde	HLD
Enzymatic Cleaner®	US	ASP® Advanced Sterilization Products	Pre-cleaner	Enzymes, surfactant	CL

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I, Aida A. Naanou, certify that the attached document which is a Self-Truthfulness Declaration for Product Registration "Acuson Sequoia" received from the company Siemens Medical Solutions, USA Inc. signed by Hyunjung Lee is an original document.

Aida A Naanou

City/County of Fairfax
Commonwealth of Virginia
Subscribed and sworn to before me
this 14th day of May, 2020
by Aida A Naanou
[Signature] Notary Public
PERLA KARA AL-SHAYEB
Reg.# 7338574 Commission expires 1/31/20 22

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- *Instruction for use for needle guides for DAX transducers*

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


GUIDANCE SYSTEM

REFERENCE GUIDE

Tracking Bracket
for use with
Siemens Healthineers
DAX
transducers

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
NOTE: Use appropriate needle length to reach target area. Verza guide channel measures 3.4cm (1.34in).

BLING TRACKING SYSTEM

With latch open, align bracket with locating features on transducer.

Close latch around transducer and click to secure. Ensure bracket is securely attached.

Place general purpose sensor onto sensor mount.

Press general purpose sensor toward bracket and click to secure.

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GUIDANCE SYSTEM

Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.


Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.

Secure cover with enclosed bands.


Inspect cover to ensure there are no holes or tears.

Remove Verza from guide storage column of gauge insert tree.

Using proper sterile technique, place Verza guide onto VerzaLink locating feature.

7. Rotate guide locking lever into locked position. Ensure guide is securely attached.

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

9. Select appropriate size gauge insert and slide into guide body chute.**BOLT-ACTION QUICK-RELEASE****1. Rotate quick-release tab to open guide insert funnel.****REMOVING GUIDANCE SYSTEM****1. Rotate guide locking lever to unlocked position to remove guide from VerzaLink locating feature of bracket.**

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REMOVING GENERAL PURPOSE SENSOR FROM BRACKET**1. Pull general purpose sensor away from bracket to remove.****CLEANING, DISINFECTING AND STERILIZING****WARNING**

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NING THE BRACKET

After each use, remove bracket from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
Inspect bracket for damage such as cracks or breakage.

INFECTING AND STERILIZING BRACKET

E: Follow manufacturer’s instructions and recommendations for concentration, time of contact and post-process procedures.

tion	Country of Origin	Manufacturer	Qualified Use	Active Ingredient	Disinfectant Type
EX® or equivalent 2% aldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Glutaraldehyde	HLD, S
EX® OPA or equivalent .55% -phthalaldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Ortho-phthalaldehyde	HLD
®	US	ASP® Advanced Sterilization Products	Pre-cleaner	Enzymes, surfactant	CL
septi® FF (new)	DE	Schülke	Soak	Succinic dialdehyde	HLD

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


NING THE GENERAL PURPOSE SENSOR

After each use remove general purpose sensor, discarding single-use, disposable components. Remove visible contaminants from general purpose sensor surface.
Inspect general purpose sensor for damage such as cracks or breakage.

INFECTING AND STERILIZING GENERAL PURPOSE SENSOR

E: Follow manufacturer’s instructions and recommendations for concentration, time of contact and post-process procedures.

tion	Country of Origin	Manufacturer	Qualified Use	Active Ingredient	Disinfectant Type
EX® or equivalent 2% aldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Glutaraldehyde	HLD, S
EX® OPA or equivalent .55% o-phthalaldehyde solutions	US	ASP® Advanced Sterilization Products	Soak	Ortho-phthalaldehyde	HLD
bl®	US	ASP® Advanced Sterilization Products	Pre-cleaner	Enzymes, surfactant	CL

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
I, Aida A. Naanou, certify that the attached document which is a Self-Truthfulness Declaration for Product Registration "Acuson Sequoia" Fusion Manual received from the company Siemens Medical Solutions, USA Inc. signed by Hyunjung Lee is an original document.

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PERLA KARA AL-SHAYEB Notary Public
Reg. # 7338574 Commission expires 1/31/20 22

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
- *Instruction for use for accessories for the imaging function in Fusion mode*

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Date

GENERAL PURPOSE SENSOR

*For use with 3D Guidance driveBAY™, driveBAY™ 2,
medSAFE™, trakSTAR™ and trakSTAR™ 2 trackers*

REFERENCE GUIDE

CE
0120

MR UNSAFE
MR-USIKKER
MR-ONVEILIG
EI MRI-TURVALLINEN
INCOMPATIBLE AVEC L'IRM
NICHT MRT-SICHER
MH ασφαλής σε MR
NON SICURO PER RM
MRで不安全

MR UNSAFE
MR-USIKKER
USO NÃO SEGURO COM RM
ИЗДЕЛИЕ НЕСОВМЕСТИМО С МРТ
核磁共振 (MR) 不安全
NO SEGURO EN RM
EJ MR-SÄKER
MR AÇISINDAN GÜVENSİZ

INDICATIONS FOR USE / INTENDED USE

This device is intended to provide physicians with tools for electromagnetic tracking of instruments with image data.

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WARNING

Prior to use inspect device for signs of damage, if damage is evident do not use.

Prior to use, ensure general purpose sensor is tracking properly.

General purpose sensor is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure general purpose sensor is properly cleaned and sterilized before each use.

To avoid possible position distortion keep general purpose sensor away from metal.

General purpose sensor is for use with CIVCO accessories only.

Refer to appropriate general purpose sensor accessory instruction for proper usage.

MR UNSAFE

Devices are MR Unsafe.

CLEANING, DISINFECTING AND STERILIZING GENERAL PURPOSE SENSOR

WARNING

Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

The recommended protocol for reprocessing the general purpose sensor has been evaluated for effectiveness and compatibility. If choosing to use a reprocessing protocol other than those listed, a validation must first be completed to ensure effectiveness and compatibility.

Do not gas sterilize or autoclave general purpose sensor.

Inspect general purpose sensor for damage such as cracks or breakage. If damage is evident, discontinue use of general purpose sensor and contact CIVCO.

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After each use remove general purpose sensor, discarding single-use, disposable components.

Remove visible contaminants from general purpose sensor surface.

Soak general purpose sensor for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL[®] Enzymatic Detergent (Johnson & Johnson). If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove general purpose sensor from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

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High-level disinfect general purpose sensor using CIDEX OPA[®] ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.

- High-level disinfect or sterilize general purpose sensor using CIDEX® Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus® (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.

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CORPORATE HEADQUARTERS:
2301 Jones Blvd., Coralville, IA 52241
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

USA OFFICE:
102 First Street South, Kalona, IA 52247-9589
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

COPYRIGHT © 2015 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

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
NEEDLE SENSOR COVER

INDICATIONS FOR USE / INTENDED USE

The cover allows use of a sensor for body surface, endocavity and intra-operative applications while helping to prevent transfer of microorganisms, body fluids and particulate material to the patient and healthcare worker.

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PRIOR TO USE

Refer to your ultrasound system's user guide for use of needle sensor.


WARNING

- Prior to use inspect device for signs of damage, if damage is evident do not use.
- Prior to use ensure needle sensor is tracking properly.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Needle sensor is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure needle sensor is properly cleaned and disinfected before each use.
- Ensure device is securely attached, failure to do so could result in inaccuracy.
- Ensure needle sensor gauge size matches eTRAX needle gauge size. Failure to do so could result in inaccuracy.
- Needle sensor system is for use with CMCO disposable needles and covers only.
- Device must be removed prior to patient defibrillation.
- To avoid possible position distortion keep needle sensor away from metal.
- eTRAX sheath is conductive and does not provide insulation. Use of ablation electrodes without an insulated cannula may result in serious injury to the patient and/or user. To ensure safe and effective use, follow the manufacturer's directions and recommended practices.
- Only stylet portion of eTRAX Needle System is a tracked instrument. After stylet is removed from sheath, the sheath should remain in pre-positioned location. If sheath is displaced at any time, reposition using tracked stylet.
- For illustration purposes only, needle sensor may be shown without a cover. Always place a cover over needle sensor to protect patients and users from cross-contamination.

SENSOR GAUGE SIZE	NEEDLE GAUGE SIZE	ACCEPTED INSTRUMENTS
12GA	12GA	≤ 14GA
14GA	14GA	≤ 16GA
16GA	16GA	≤ 18GA
18GA	18GA	≤ 20GA

BLING NEEDLE SYSTEM

Using proper sterile technique, insert needle sensor through opening in cover, taking care to minimize bending of needle sensor.

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

NG

Prior to use ensure sheath is securely attached to stylet hub, do not puncture without sheath in place.

Using proper sterile technique, insert needle sensor into needle, taking care to minimize bending of needle sensor.

Lock needle onto needle handle, taking care to avoid puncturing cover.


T REMOVAL

Once target is reached, rotate stylet hub out of sheath hub by turning counter-clockwise.

CLEANING AND DISINFECTING NEEDLE SENSOR**WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
 - *The recommended protocol for reprocessing the needle sensor has been evaluated for effectiveness and compatibility. If choosing to use a reprocessing protocol other than those listed, a validation must first be completed to ensure effectiveness and compatibility.*
 - *Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.*
 - *Ensure needle sensor does not bend during cleaning and disinfecting.*
 - *Do not gas sterilize or autoclave needle sensor.*
 - *Inspect needle sensor for damage such as cracks or breakage. If damage is evident, discontinue use of needle sensor and contact CIVCO.*
1. After each use, remove needle sensor from needle, discarding single-use, disposable components. Remove visible contaminants from needle sensor surface.
 2. Clean and disinfect all surfaces with common germicidal or antiseptic such as Sani-cloth® HB, CIDEX OPA® (Johnson & Johnson), ENZOL® (Johnson & Johnson), T-Spray™, Protex™ or alcohol.
 3. Follow manufacturer's directions for all cleaning and disinfection.

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MPS GmbH
Borngasse 20
35619 Braunfels, Germany
Fax: 49 6442 32578

Australian
Sponsor Address

Emergo Australia
Level 20
Tower II, Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

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COPYRIGHT © 2017 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ETRAK IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

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COPYRIGHT © 2016 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. VERZA AND VERZALINK ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
COPYRIGHT © 2018 ALL RIGHTS RESERVED. CNCO IS A REGISTERED TRADEMARK OF CNCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.
COPYRIGHT © 2019 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
COPYRIGHT © 2018 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. VERZA AND VERZALINK ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.