Medtronic

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
	Хаб	30010, 30012, 30014, 30011, 30055, 30050	Полипропил ен (голубой)	100	Кратковременный (< 24 часов), с циркулирующей кровью																																												
		30110, 30114, 30112, 30212, 30155, 30255	Полистирен (голубой)	100	Кратковременный (< 24 часов), с циркулирующей кровью																																												
	<table><tr><th colspan="2">Component</th><th>Component material</th><th>Percentage ratio (%)</th><th>Type of contact with patient/person</th></tr><tr><td rowspan="3">Basket tip 30010, 30012, 30014</td><td>End of the coronary ostial perfusion cannula</td><td>PVC - blue with BaSO₄</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td>Glue</td><td>Cyclohexanone</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td>Coupling</td><td>ABS Plastic (blue)</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td rowspan="3">Spherical tip 30011</td><td>End of the coronary ostial perfusion cannula</td><td>PVC - blue with BaSO₄</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td>Glue</td><td>Cyclohexanone</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td>Coupling</td><td>ABS Plastic (blue)</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td rowspan="3">Soft convex tip 30055</td><td>Tube - Shrink</td><td>Polyester</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td>Part A</td><td>Silicone</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr><tr><td>Part B</td><td>Silicone</td><td>100</td><td>Short-term (<24 hrs) with circulating blood</td></tr></table>					Component		Component material	Percentage ratio (%)	Type of contact with patient/person	Basket tip 30010, 30012, 30014	End of the coronary ostial perfusion cannula	PVC - blue with BaSO ₄	100	Short-term (<24 hrs) with circulating blood	Glue	Cyclohexanone	100	Short-term (<24 hrs) with circulating blood	Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood	Spherical tip 30011	End of the coronary ostial perfusion cannula	PVC - blue with BaSO ₄	100	Short-term (<24 hrs) with circulating blood	Glue	Cyclohexanone	100	Short-term (<24 hrs) with circulating blood	Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood	Soft convex tip 30055	Tube - Shrink	Polyester	100	Short-term (<24 hrs) with circulating blood	Part A	Silicone	100	Short-term (<24 hrs) with circulating blood	Part B	Silicone	100	Short-term (<24 hrs) with circulating blood
	Component		Component material	Percentage ratio (%)	Type of contact with patient/person																																												
	Basket tip 30010, 30012, 30014	End of the coronary ostial perfusion cannula	PVC - blue with BaSO ₄	100	Short-term (<24 hrs) with circulating blood																																												
		Glue	Cyclohexanone	100	Short-term (<24 hrs) with circulating blood																																												
		Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood																																												
	Spherical tip 30011	End of the coronary ostial perfusion cannula	PVC - blue with BaSO ₄	100	Short-term (<24 hrs) with circulating blood																																												
		Glue	Cyclohexanone	100	Short-term (<24 hrs) with circulating blood																																												
Coupling		ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood																																													
Soft convex tip 30055	Tube - Shrink	Polyester	100	Short-term (<24 hrs) with circulating blood																																													
	Part A	Silicone	100	Short-term (<24 hrs) with circulating blood																																													
	Part B	Silicone	100	Short-term (<24 hrs) with circulating blood																																													

Component		Component material	Percentage ratio (%)	Type of contact with patient/person
Basket tip 30010, 30012, 30014	End of the coronary ostial perfusion cannula	PVC - blue with BaSO ₄	100	Short-term (<24 hrs) with circulating blood
	Glue	Cyclohexanone	100	Short-term (<24 hrs) with circulating blood
	Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
Spherical tip 30011	End of the coronary ostial perfusion cannula	PVC - blue with BaSO ₄	100	Short-term (<24 hrs) with circulating blood
	Glue	Cyclohexanone	100	Short-term (<24 hrs) with circulating blood
	Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
Soft convex tip 30055	Tube - Shrink	Polyester	100	Short-term (<24 hrs) with circulating blood
	Part A	Silicone	100	Short-term (<24 hrs) with circulating blood
	Part B	Silicone	100	Short-term (<24 hrs) with circulating blood

	Soft concave tip 30050	Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
		Tube - Shrink	Polyester	100	Short-term (<24 hrs) with circulating blood
		Part A	Silicone	100	Short-term (<24 hrs) with circulating blood
		Part B	Silicone	100	Short-term (<24 hrs) with circulating blood
		Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
	90° angle basket tip 30110, 30112, 30114	Coronary Bending 90°	Acryl (white)	100	Short-term (<24 hrs) with circulating blood
		End of the coronary ostial perfusion cannula	PVC - blue with BaSO4	100	Short-term (<24 hrs) with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
	45° angle basket tip 30212	End of the coronary ostial perfusion cannula	PVC - blue with BaSO4	100	Short-term (<24 hrs) contact with circulating blood
		Coronary Bending 45°	Acryl white	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
	90° angle soft tip 30155	Primer	Silicone	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Transparent in-line sealer	100	Short-term (<24 hrs) contact with circulating blood
		End of the coronary ostial perfusion cannula - soft	Silicone	100	Short-term (<24 hrs) with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
		Connector 90°	Acryl white	100	Short-term (<24 hrs) with circulating blood

	45° angle soft tip 30255	Primer	Silicone	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Transparent in-line sealer	100	Short-term (<24 hrs) contact with circulating blood
		End of the coronary ostial perfusion cannula - soft	Silicone	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
		Connector 45°	Acryl white	100	Short-term (<24 hrs) contact with circulating blood
	Soft bulb beveled tip 30315, 30317, 30320		Silicone, gel	100	Short-term (<24 hrs) with circulating blood
	Silicone body 30315, 30317, 30320		Silicone, gel	100	Short-term (<24 hrs) with circulating blood
	Malleable stainless steel tube 30010, 30012, 30014, 30011, 30055, 30050, 30110, 30114, 30112 30212, 30155, 30255		Stainless steel	100	Short-term (<24 hrs) with circulating blood
	Female luer lock connector	30010, 30012, 30014, 30011, 30055 30050	Polypropylene (blue)	100	Short-term (<24 hrs) with circulating blood
		30110, 30114, 30112, 30212, 30155, 30255	Acryl white	100	Short-term (<24 hrs) with circulating blood
		30315, 30317, 30320	Polypropylene	100	Short-term (<24 hrs) with circulating blood
	Plastic handle 30110, 30114, 30112, 30212, 30155, 30255		Acryl white	100	Short-term (<24 hrs) with circulating blood

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Information in Russian sticker

The following information should be applied on a Russian-language sticker (word order and text font may vary)

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The individual packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Sign of conformity in Russia;
- Warning and informing symbols and inscriptions (if applicable).

Sample of the sticker

 Non-toxic   EC REP	Name of medical device	
	Manufacturer: Medtronic Inc.	Code UPN
	Country of origin: USA	Code CFN
	Manufacturer's authorized representative: LLC «Medtronic»	
	Registration Certificate № _____ or _____	

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Technical parameters:

1. Outer diameter of the end of the clear body with one-way valve and 3-mm beveled tip: 3 mm

2. Cannula total length: 6.1 cm


3. Cannula effective length: 38 mm

4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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Technical parameters:

1. Outer diameter of the end of the clear body with one-way valve and 3-mm blunt tip: 3 mm

2. Cannula total length: 6 cm

3. Cannula effective length: 38 mm


4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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Technical parameters:

1. Outer diameter of the end of the vessel cannula with clear body and 3-mm blunt tip: 3 mm

2. Cannula total length: 5.5 cm


3. Cannula effective length: 31 mm

4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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Technical parameters:

1. Outer diameter of the end of the vessel cannula with clear body and 2-mm blunt tip: 2 mm

2. Cannula total length: 5.5 cm

3. Cannula effective length: 32 mm


4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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Technical parameters:

1. Outer diameter of the end of the vessel cannula with clear body and 4-mm acorn tip: 4 mm

2. Cannula total length: 5.5 cm


3. Cannula effective length: 33 mm

4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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Technical parameters:

1. Outer diameter of the end of the vessel cannula with radiopaque body and 3-mm beveled tip: 3 mm

2. Cannula total length: 5.4 cm


3. Cannula effective length: 31 mm

4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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Technical parameters:

1. Outer diameter of the end of the clear body with 3-mm beveled tip: 3 mm


2. Cannula total length: 6.1 cm

3. Cannula effective length: 38 mm


4. Connection dimension of the locking female luer fitting:

Inner diameter (locking female luer fitting): 4,315 (+0 / -0,005) mm

Outer diameter (locking female luer fitting): 7,83 (+0 / -0,1) mm

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
	Radiopaque body with one-way valve and 3-mm blunt tip 30021	One - way check valve - "duck nose" type	Silicone	100	Short-term (<24 hrs) with circulating blood
		Vessel cannula body - radiopaque	PVC - white with BaSO4	100	Short-term (<24 hrs) with circulating blood
	Clear body with radiopaque band and 3-mm blunt tip 30022	Glue	Cyclohexanone, PVC	100	Short-term (<24 hrs) contact with intact skin
		Expander-downward arc	Tecoflex, BaSO4 40%	100	Short-term (<24 hrs) with circulating blood
		Clear body with 3-mm blunt tip	PVC transparent	100	Short-term (<24 hrs) with circulating blood
	Radiopaque body with 3-mm blunt tip 30023		PVC - white with BaSO4	100	Short-term (<24 hrs) with circulating blood
9. Данные маркировке медицинского изделия и его упаковке/ Labelling and Packaging of medical device	<u>Маркировка индивидуальной упаковки: Individual package labeling:</u> <u>Для канюль для перфузии устья коронарной артерии: / For coronary ostial perfusion cannulas:</u> 				

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

The final packaging for the Russian market is accompanied with a sticker.

Devices supplied in packaging marked with the following information:

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	Name of medical device	
Non-toxic 	Manufacturer: Medtronic Inc.	Code UPN
	Country of origin: USA	Code CFN
	Manufacturer's authorized representative: LLC «Medtronic»	
	Registration Certificate № ____ от ____	
	Date of manuf../Lot/UBD: See on pack	See IFU See packaging symbols

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
1. Retractor length: 45 mm ((+0.01 mm (+-0.0004 inch))
2. Retractor width: 29 mm ((+0.01 mm (+-0.0004 inch))
3. The inner diameter of the tube is 0.635 mm; 0.025 inch (+- 0.05 mm (± 0.002 inch))
4. The outer diameter is 1.2 mm; 0.043 inches (+- 0.05 mm (±0.002 inches))

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Retractor for coronary artery exposure 16130, 16150

Component	Component material	Percentage ratio (%)	Type of contact with patient/person
Stainless steel wire	Stainless steel	100	Short-term (<24 hrs) with circulating blood
Soft, polymer coating	extruded radiopaque polyester-polyurethane elastomer	100	Short-term (<24 hrs) with circulating blood
Teeth (available in 3.0 and 5.0 mm length)	Stainless steel	100	Short-term (<24 hrs) with circulating blood
Substrate	polypropylene synthetic paper	100	Short-term (<24 hrs) with circulating blood

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

Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby swear (or affirm) that the attached documentation is true and unaltered copies of the original document.

Name: Linda Julia Davilla
Title: Regulatory Affairs Specialist

____ May 29, 2018 ____
Date

~~Haus~~
May 29, 2018

DLP®

Vessel Cannula
Canule
Gefäßkanüle
Cánula vascular
Vaatskanule
Cannula vascolare
Venekanyla
Verisuonikanyyli
Kärlkanyl
Κάνουλα αγγείων
Venekanyla
Kaniula naczyniowa
Cánula de vasos
Damar Kanülü
Cévní kanyla
Vénakanül
Канюля для сосудов

How
May 29, 2018

Instruccions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
Bruksanvisning • Οδηγίες χρήσης • Brugsanvisning • Instrukcja użytkowania
Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati útmutató
Инструкция по эксплуатации

USA Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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Vessel Cannula

Product Description

These vessel cannulas have soft tips and locking female luer. Several models feature a one-way valve that prevents backflow through the cannula. Some models offer a radiopaque body or band (see Figure 1). Sterile, nonpyrogenic, single use.

Indications for Use

These cannulas are intended for use up to six hours or less in conjunction with cardiopulmonary bypass surgery or in vascular surgery to perfuse a vein graft or to help check for leaks in a harvested vein which will be used for a graft.

Contraindications

This device is not intended for use except as indicated above.

Warnings

Extreme caution is needed when using a vessel cannula to check for leaks in a harvested vein. The vein could rupture under conditions of extreme pressure.

Precautions

Note: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, and the type of surgical procedure being performed. When infusing cardioplegia solution, use the lowest possible pressures consistent with surgical techniques to minimize hemolysis, or damage to the vessels.

Purge all the air from the vessel cannula before connection.

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events, and dislocation. Vessel damage and complications at the puncture site may occur if the instructions for Use are not followed.

Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique.
2. If using the vessel cannula to check for leaks in a harvested vein, connect the cannula to a syringe filled with sterile saline solution. Tie off one end of the vein and insert the tip of the vessel cannula into the open end of the vein. Tie a piece of suture around the vein and the cannula to keep the cannula in the vein. Gently fill the vein with sterile saline solution and check for any leaks.
3. If using the vessel cannula to perfuse an artery, connect the vessel cannula to the administration set and prime the cannula to purge any air that may be in the cannula. Place the cannula in the artery or proximal end of the vein graft and flow at the lowest possible pressure that will adequately perfuse the artery.
4. At the completion of the procedure, disconnect the cannula and dispose of it according to hospital policy.

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Disclaimer of Warranty

ALTHOUGH THE VESSEL CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Disclaimer of Warranty

ALTHOUGH THE VESSEL CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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2006/OCT/02 at 2:01 p.m. Doc number: M835132A001 Rev. A (NO17)

Medtronic

Advancing Life • Restoring Health • Sustaining Life

Europe

Europe/Africa/Middle East Headquarters
Medtronic International Trading Sarl
Route du Mollet 31
Case Postale
CH - 1131 Tolochenaz
Switzerland
Internet: www.medtronic.co.uk
Tel. 41-21-802-7000
Fax 41-21-802-7900

Medtronic E.C. Authorized
Representative/Distributor by
Medtronic B.V.
Earl Bakkenstraat 10
8432 PJ Heerlen
The Netherlands
Tel. 31-45-606-8000
Fax 31-45-606-8008

Asia-Pacific

Japan
Medtronic Japan
Comodio Shiodome SF
2-14-1 Higashi-Shimbashi, Minato-ku
Tokyo 105-0021
Japan
Tel. 81-3-6430-2011
Fax 81-3-6430-7140

Australia
Medtronic Australasia Pty. Ltd.
Business Address:
87 Waterloo Road
North Ryde NSW 2113
Mailing Address:
PO Box 945
North Ryde NSW 1670
Australia
Tel. 61-2-9857-8000
Fax 61-2-9878-5100

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Arteriotomy Cannula

Product Description

The IMA (Internal Mammary Artery) Cannula consists of a plastic female slip fit connector and 25-gauge stainless steel tubing that ends in a 1-cm bulb. The Arteriotomy Cannulae consist of a plastic winged female luer lock connector with a length of polyurethane tubing that ends in a radiopaque bulb tip.

Sterile, nonpyrogenic, single use.

Indications for Use

These cannulae are intended for the perfusion of cardioplegic solution directly into a cardiac bypass graft or artery, or retrograde flushing of the free internal mammary artery with a vasodilator. These devices are intended for use up to six hours or less.

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Warnings

These cannulae may not be detectable on x-ray in some cases due to the variances in anatomical densities and x-ray techniques.

Precautions

Notes: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, and the type of surgical procedure.

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Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events, and dislocation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique.
2. Connect the cannula to the cardioplegic solution delivery tubing or a pressure sensing syringe.
3. Prime the cannula to remove any air before perfusing the bypass graft or artery.
4. Carefully insert the cannula into the target vessel using caution to avoid any trauma to the internal intima.
5. Slowly increase the flow to the desired rate or the maximum line pressure prescribed by the surgeon. When using a pressure sensing syringe, caution must be exercised so as not to over pressurize or over inflate the artery or graft.
6. When the surgical procedure is completed, dispose of the cannula according to hospital policy. **DO NOT REUSE.**

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Disclaimer of Warranty

ALTHOUGH THE ARTERIOTOMY CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining

portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Disclaimer of Warranty

ALTHOUGH THE ARTERIOTOMY CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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2009/APR/21 at 4:34 p.m. Doc number: MS34522A001 Rev. A [NO17]

Medtronic

Measuring Data. Powering Health. Saving Lives.

Europe

Europe/Africa/Middle East Headquarters

Medtronic International Trading Ltd

Rue du Mollet 31

Case Postale

CH - 1131 Tolochenaz

Switzerland

Internet: www.medtronic.co.uk

Tel. 41-21-802-7000

Fax 41-21-802-7800

Medtronic E.C. Authorized

Representative/Distributor by

Medtronic B.V.

Earl Bakkenstraat 18

6422 PJ Heerlen

The Netherlands

Tel. 31-45-666-8000

Fax 31-45-666-8000

Asia-Pacific

Japan

Medtronic Japan

Comedia Shinjuku 5F

2-14-1 Higashi-Shinjuku, Minato-ku

Tokyo 105-8021

Japan

Tel. 81-3-6430-2011

Fax 81-3-6430-7140

Australia

Medtronic Australia Pty. Ltd.

Business Address:

97 Waterloo Road

North Ryde NSW 2113

Mailing Address:

PO Box 946

North Ryde NSW 1670

Australia

Tel. 61-2-9867-6000

Fax 61-2-9878-6100

Asia

Medtronic International Ltd.

Suite 1602 16/F, Manilla Plaza

The Lee Gardens, 33 Hysan Avenue

Causeway Bay

Hong Kong

Tel. 852-2981-4088

Fax 852-2981-0313

Americas

Latin America

Medtronic Latin America

3750 NW 87th Avenue

Suite 700

Miami, FL 33178

USA

Tel. 305-600-6328

Fax 788-708-4244

Canada

Medtronic of Canada Ltd.

6733 Kilmer Road

Markham, Ontario L3R 1W3

Canada

Tel. 905-826-8020

Fax 905-826-8820

Toll-free in Canada:

1-800-388-6346

United States

Manufacturer:

Medtronic, Inc.

710 Medtronic Parkway

Minneapolis, MN 55432

USA

Internet: www.medtronic.com

Tel. 763-514-4000

Fax 763-381-8100

Toll-free in the USA:

1-800-328-2518

(24-hour consultation service)

© Medtronic, Inc. 2008

All Rights Reserved

M834622A001 Rev. A

DLP®

Coronary Ostial Perfusion Cannula
Canule de perfusion coronaire ostiale
Perfusionskanüle für Koronarostien
Cánula de perfusión para ostium coronario
Perfusiecanule voor coronaire ostie
Cannula per perfusione coronarica ostiale
Kanyla for koronar-ostial perfusion
Syöttimen ostiaalinen perfusiokanyyli
Kanyl för koronar-ostial perfusion
Κάνουλα οστέιας στεφανιαίας αρτηρίας
Perfusionskanyla til det koronare ostium
Kaniula do kardiolegii dowieńcowej
Cânula de perfusão do ostium coronário
Koroner Ostial Perfúzyon Kanüüü
Ostiaální perfúzní kanyla pro koronární intervenci
Koronáriszűjtődékba illeszthető perfúziós kanül
Канюля для перфузии устья коронарной артерии

How
May 29, 2018

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
Bruksanvisning • Οδηγός χρήσης • Brugsanvisning • Instrukcja użytkowania
Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati útmutató
Инструкция по применению

(USA) Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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DLP®

Coronary Ostial Perfusion Cannula

Product Description

These coronary cannulae consist of a basket-style tip, soft (silicone) tip, or spherical tip which is attached to a malleable stainless steel tube, or a soft bulb beveled tip with an integral silicone body. The cannulae terminate with a locking female luer fitting. The French diameter is measured at the base of the tip adjacent to the sealing flange or across the diameter of the soft bulb beveled tip style.

Sterile, nonpyrogenic, single use.

Indications for Use

These cannulae are intended for use in conjunction with cardiopulmonary bypass surgery up to six hours or less for delivery of cardioplegia solutions directly to the coronary arteries.

Contraindications

This device is not intended for use except as indicated above. If the ostium is calcified, the basket tip cannula, the soft bulb beveled tip cannula, and spherical tip cannula models should not be used.

Warnings

Utilization of this cannula may not result in adequate perfusion of the heart if the coronary arteries have extensive disease and/or occlusions.

Infusion pressures **must be carefully monitored and controlled** so they do not exceed the maximum level specified by the surgeon. If the pressure level is exceeded, the coronary arteries may be damaged.

Do not use a basket tip cannula, a soft bulb beveled tip cannula, or a spherical tip cannula that is larger in size than the coronary ostium. Compression of the tip could restrict flow.

Do not insert the flange of the basket tip cannula, the bulb of the beveled tip cannula, or the large diameter bulb of the spherical tip cannula into the ostium as this may block the lateral branches near the ostium resulting in inadequate perfusion of these branches.

Precautions

Note: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, and the type of surgical procedure being performed.

When infusing cardioplegia solution use the lowest possible pressures consistent with the surgical techniques to minimize hemolysis, or damage to the vessels.

These cannulae are designed and intended for single use only. **DO NOT REUSE.**

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events and dislocation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique.
2. Select the correct size and style of cannula. Connect the cannula to a cardioplegia delivery system which has an integral means of monitoring the line pressure.
3. Prime the cannula and flush to remove any air before cannulation.
4. The cannula shaft may be bent into the desired shape. Avoid any sharp kinks.
5. With the cardioplegia solution flowing at a very low rate, fill cup of cannula outlet (basket tip model only) and place over the coronary ostium being careful not to entrap air. With the soft bulb beveled tip cannula, forceps may be applied proximal to the bulb to aid in placement/retention of the cannula tip into the coronary ostium. Increase the flow rate to the desired rate or the maximum line pressure prescribed by the surgeon. Infuse the prescribed amount of cardioplegia solution (Figure 2).
6. Remove the cannula.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE CORONARY OSTIAL PERFUSION CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS,

MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE CORONARY OSTIAL PERFUSION CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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

Medtronic

Operational Head Quarters
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com

tel 800.633.8766

Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby swear (or affirm) that the attached documentation are true and unaltered copies of the original document.

Diane Howell


Senior Regulatory Affairs Specialist | Cardiac Surgery, CSH

Medtronic

8200 Coral Sea Street NE, MVS 83 | Mounds View, MN, 55112 | USA

1-September 2020

Date

As a registered Notary Public for the State of Minnesota,

I hereby certify that this is a true and accurate copy of the

original document. The original is filed at Medtronic, Inc.

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			Tetrahydrofuran		
		Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
	Soft convex tip 30055	Tube - Shrink	Polyester	100	Short-term (<24 hrs) with circulating blood
		Part A	Silicone	100	Short-term (<24 hrs) with circulating blood
		Part B	Silicone	100	Short-term (<24 hrs) with circulating blood
		Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
		Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
	Soft concave tip 30050	Tube - Shrink	Polyester	100	Short-term (<24 hrs) with circulating blood
		Part A	Silicone	100	Short-term (<24 hrs) with circulating blood
		Part B	Silicone	100	Short-term (<24 hrs) with circulating blood
		Coupling	ABS Plastic (blue)	100	Short-term (<24 hrs) with circulating blood
	90° angle basket tip 30110, 30112, 30114	Coronary Bending 90°	Acryl (white)	100	Short-term (<24 hrs) with circulating blood
		End of the coronary ostial perfusion cannula	PVC - blue with BaSO4	100	Short-term (<24 hrs) with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
	45° angle basket tip 30212	End of the coronary ostial perfusion cannula	PVC - blue with BaSO4	100	Short-term (<24 hrs) contact with circulating blood
		Coronary Bending 45°	Acryl white	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood

	90° angle soft tip 30155	Primer	Silicone	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Transparent in-line sealer	100	Short-term (<24 hrs) contact with circulating blood
		End of the coronary ostial perfusion cannula - soft	Silicone	100	Short-term (<24 hrs) with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
		Connector 90°	Acryl white	100	Short-term (<24 hrs) with circulating blood
	45° angle soft tip 30255	Primer	Silicone	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Transparent in-line sealer	100	Short-term (<24 hrs) contact with circulating blood
		End of the coronary ostial perfusion cannula - soft	Silicone	100	Short-term (<24 hrs) contact with circulating blood
		Adhesive	Cyanoacrylate adhesive	100	Short-term (<24 hrs) contact with circulating blood
		Connector 45°	Acryl white	100	Short-term (<24 hrs) contact with circulating blood
	Soft bulb beveled tip 30315, 30317, 30320		Silicone, gel	100	Short-term (<24 hrs) with circulating blood
	Silicone body 30315, 30317, 30320		Silicone, gel	100	Short-term (<24 hrs) with circulating blood
	Malleable stainless steel tube 30010, 30012, 30014, 30011, 30055, 30050, 30110, 30114, 30112 30212, 30155, 30255		Stainless steel	100	Short-term (<24 hrs) with circulating blood

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