«УТВЕРЖДАЮ» / “APPROVE”
ETHICON ENDO-SURGERY, LLC, USA

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**INSTRUCTION FOR USE
OF THE MEDICAL DEVICE**

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
2018 year / г.

[Signature]
27 Feb 2018
SARAH WARD
NOTARY PUBLIC
STATE OF OHIO
Comm. Expires
12-02-2022

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Important: This package insert is designed to provide instructions for use of the ECHELON™Flex Articulating Endoscopic Linear Cutters. It is not a reference to surgical techniques.



ENDOPATH and ECHELON are trademarks of Ethicon Endo-Surgery.

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
portion of the shaft to pivot to facilitate lateral access to the operative site.

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Caution: Do not load the instrument more than 12 times for a maximum of 12 firings per instrument. The use of the instrument with staple line reinforcement material may reduce the number of firings.

ECHELON Flex Instrument Product Codes

Instrument Code	Description	Shaft Length
EC45A	ECHELON Flex 45 Articulating Endoscopic Linear Cutter	340 mm
SC45A	ECHELON Flex 45 Compact Articulating Endoscopic Linear Cutter	280 mm
EC45AL	ECHELON Flex 45 Long Articulating Endoscopic Linear Cutter	440 mm
EC60A	ECHELON Flex 60 Articulating Endoscopic Linear Cutter	340 mm
SC60A	ECHELON Flex 60 Compact Articulating Endoscopic Linear Cutter	280 mm
LONG60A	ECHELON Flex 60 Long Articulating Endoscopic Linear Cutter	440 mm

ECHELON Reload Product Codes

Note: The ECHELON Flex instruments may only be used with the ECHELON reloads.

Reload code	Closed Staple Height	Staple Line Length	Reload Color	Number of Staplers	Stapler Rows
ECR45M	0,75 mm	45 mm	Gray	70	6
ECR45W	1,0 mm	45 mm	White	70	6
ECR45B	1,5 mm	45 mm	Blue	70	6

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Clinical Section

Examples of procedures in which linear cutters and staplers may be used include:

General	Гинекология	Урология	Торакальная хирургия
- Anastomoses - Partial Gastroplasty - Appendectomy	Vaginal Hysterectomy	- Prostatectomy - Nephrectomy - Cystectomy	- Lung Resection and Biopsy - Pulmonary Wedge Resection

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Instructions for Use

Prior to Using Instrument:

- 1 Verify that the reload size matches the instrument size to be used (e.g. use an ECHELON 45 reload with an ECHELON Flex 45 instrument).
- 2 Become familiar with the articulating action of the instrument by pulling on the fins of the rotating knob with the index finger and articulating the jaws with the other hand. (Illustration 2) As lateral pressure is applied in either direction, the jaws pivot to a maximum angle of 45 degrees, at which point the pivoting force increases to signal the end of the articulating arc. Maintain the lateral pressure while releasing the rotating knob to lock the angle of the shaft. After locking the jaws, release the lateral pressure. The jaws will also lock at approximate 15 and 30 degrees by releasing the rotating knob before the end of the articulation arc. To return the jaws to the straight position, pull on the fins of the rotating knob to release the lock; the jaws will return automatically.
- 3 Verify compatibility of all instruments and accessories prior to using the instrument (refer to **Warnings and Precautions**).

Loading the Instrument:

- 4 Using sterile technique, remove the instrument and reload from their respective packages. To avoid damage, do not flip the instrument or reload into the sterile field.
- 5 Remove the multiple stroke firing instructions tag and plastic sleeve from the instrument. To do this, slide the tag and sleeve over the shaft of the instrument.
- 6 Prior to loading, ensure the instrument is in the open position. (Illustration 3)
- 7 Examine the reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

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Caution: Tissue thickness should be carefully evaluated prior to using the instrument.

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- 8 Insert the new reload by sliding it against the bottom of the cartridge jaw until the cartridge alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now loaded and ready for use. (Illustration 4)

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Using the Instrument:

- 9 Close the jaws of the instrument by squeezing the closing trigger (light gray and labeled with a Number "1") until it locks in place. An audible click indicates that the closing trigger is locked. (Illustration 5)

Caution: Do not squeeze the firing trigger (dark gray and labeled with a Number "2") at this time.

- 10 Introduce the instrument into the body cavity through a trocar of the appropriate size or through an incision. (Illustration 6) When using a trocar, the instrument jaws must be visible past the trocar sleeve before opening the jaws.

Caution: When placing the instrument through the trocar or incision, avoid advancing the firing trigger (2) accidentally. If this occurs, the instrument will lockout and the stroke indicator will display a lock symbol; in this state, the instrument will not allow the

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anvil release button to reopen the instrument jaws. To recover from this condition, remove the instrument, push the red manual knife reverse switch downward to reverse the knife motion; the knife indicator will display an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. Squeeze the firing trigger (2) completely until it rests on the closing trigger (1); the stroke count indicator will display zero to indicate the knife has been returned to its home position. (Illustration 7) Press the anvil release button and replace the reload, then follow the instructions above for Loading the Instrument.

- 11 Once in the cavity, press the anvil release button to reopen the instrument jaws and return the closing trigger (1) to its original position. (Illustration 8)
- 12 If needed, rotate the jaws by pushing on the fins of the rotating knob with the index finger using a downward or upward pressure. The instrument shaft will rotate freely in either direction.
- 13 To articulate the jaws inside the body cavity, select an appropriate surface (a body structure, organ, or another instrument) to use as a grounding surface to bend the jaws, ensuring that it stays within the field of view. Rest the lateral side of the jaws that are opposite to the desired direction of articulation against the grounding surface (the jaws must be open in order to articulate the instrument). (Illustration 9)

WARNING: Do not attempt to articulate by pressing the front of the jaws against the grounding surface as tissue damage or tissue trauma may occur. (Illustration 10)

Pull the fins of the rotating knob back with the index finger and apply a sweeping motion towards the side articulation is wanted while gently pushing the instrument handle towards the grounding surface. Maintain the jaws pressed against the grounding surface during this action. Once the desired articulation angle is reached, release the rotating knob to lock the angle (the instrument will only lock at pre-determined angles – 15°, 30° and 45°).

Caution: The instrument can only achieve a maximum articulation angle of 45°. When using body structures or organs as a grounding surface, particular attention should be placed to the visual cues and tactile feedback received from the instrument. When the maximum angle is reached, the force will increase indicating the maximum

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angle has been reached. Avoid applying excessive pressure to the tissue as tissue damage or tissue trauma may occur.

14 Position the instrument around the tissue to be stapled.

Caution: Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line.

The distal black lines on the anvil and cartridge channel designate the ends of the staple line. The line on the cartridge channel that reads “cut” references the cut line on the device. (Illustration 1)

Caution: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.

15 After positioning the instrument jaws, close the jaws by squeezing the closing trigger (1) toward the handle until it locks. An audible/tactile click indicates that the closing trigger and the jaws are locked. (Illustration 5)

Caution: Ensure tissue has not extended (extruded) proximal to the proximal black line on the instrument. Tissue forced into the instrument proximal to the black line may be transected without staples.

Caution: If the closing trigger (1) is difficult to lock, **reposition the instrument and take a smaller amount of tissue.**

Ensure that the proper reload selection has been made. (Refer to the **ECHELON Reload Product Codes table.**)

Caution: If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use instrument.

Note: When firing across thick tissue, holding the jaws in place for 15 seconds after closing and prior to firing may result in better compression and staple formation.

16 Fire the instrument. The complete firing action requires 4 complete strokes to cut the full reload length and return the knife to its home position. For each stroke, squeeze and release the firing trigger (2) completely with a continuous, smooth stroke until it rests on the closing trigger (1). (Illustration 11) The numbers on the stroke count

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indicator will advance with each stroke. After the third stroke, the knife direction indicator will display an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. At the end of the fourth stroke, the stroke count indicator will display "0" to indicate the knife has returned to its home position. The status of the transection can also be determined by observing the knife blade indicator throughout the firing action.

Note: Once the firing cycle has been initiated, it should be completed by completing all four strokes of the firing cycle to return the knife to the home position. If it is necessary to interrupt the firing sequence, push the red manual knife reverse switch downward to reverse the knife motion; the knife direction indicator will display an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. To complete, squeeze the firing trigger (2) completely until it rests on the closing trigger (1), after which the stroke count indicator will display "0" to indicate the knife has returned to its home position.

Caution: Attempting to force the trigger to complete the firing stroke with too much tissue between the jaws, or with dense/ thick tissue between the jaws, may result in increased forced to fire or instrument failure.

Caution: If the instrument locks out, the stroke count indicator will display a lockout symbol. Stop and push the red manual knife reverse switch downward to reverse the knife motion; the knife direction indicator will display an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. Squeeze firing trigger (2) completely until it rests on the closing trigger (1). The stroke count indicator will display "0" to indicate the knife has returned to its home position. Then press the anvil release button to release the jaws from the tissue. Close the instrument jaws by pressing the closing trigger (1), remove the instrument, and replace the cartridge. Firing through the lockout mechanism will break the instrument.

Caution: The use of staple line reinforcement material with the instrument may require an increased force to fire and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.

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17 Press the anvil release button to separate the instrument jaws and allow both triggers to return to their original positions. (Illustration 8)

Caution: If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the retracted position by verifying that the stroke count indicator displays "0" and knife direction indicator points towards the proximal side of the instrument, or that the knife blade indicator is in the home position. If the stroke count indicator or knife blade indicator is not in the home position, push the red manual knife reverse switch downward to reverse the knife motion and squeeze the firing trigger (2) completely until it rests on the closing trigger (1). Press the anvil release button. If the jaws do not open at this point, then gently pull the closing trigger (1) upward (away from the handle) until both firing and closing triggers return to their original positions.

Caution: If the firing mechanism becomes inoperative, **do not** continue to use the instrument.

18 Gently pull the instrument away from the transected tissue and ensure it is released from the jaws. Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with electrocautery, manual sutures or other appropriate techniques.

19 Before removing articulating instruments, move the jaws away from any obstruction inside the body cavity while keeping the jaws open and within the field of view, and pull on the fins of the rotating knob. The jaws will return to the straight position automatically.

20 To remove the instrument from the cavity, squeeze the closing trigger (1) until it locks, closing the jaws. Completely withdraw the instrument in the closed position. (Illustration 12)

Reloading the Instrument:

21 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

22 Prior to reloading, ensure the instrument is in the open position. (Illustration 3)

23 Push upward (toward the anvil) to unsnap the reload from the cartridge jaw. Discard the used reload. (Illustration 13)

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WARNING: Prior to reloading the instrument, rinse the anvil and cartridge jaw in sterile solution and then wipe the anvil and cartridge jaw to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and cartridge jaw.

- 24 Examine the new reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

Caution: Selection of the appropriate staple cartridge should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to fire and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.

WARNING: Loading the incorrect reload for the instrument size or model may result in transecting the tissue without sealing (e.g. loading a 60 mm reload into a 45 mm instrument).

- 25 Insert the new reload by sliding it against the bottom of the cartridge jaw until the cartridge alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use. (Illustration 14)

Caution: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)

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Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- When using other technologies (e.g., electrocautery) in the procedure, observe the precautions suggested by the original equipment manufacturer to avoid the hazards associated with their use.
- Loading the incorrect reload for the instrument size or model may result in transecting the tissue without sealing (e.g. loading a 60 mm reload into a 45 mm instrument).
- Do not attempt to articulate by pressing the front of the jaws against the grounding surface as tissue damage or tissue trauma may occur.
- Prior to reloading the instrument, rinse the anvil and cartridge jaw in sterile solution and then wipe the anvil and cartridge jaw to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and cartridge jaw.
- Before using, remove the staple retaining cap and observe the surface of the reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)
- When placing the instrument through the trocar or incision, avoid advancing the firing trigger (2) accidentally. If this occurs, the instrument will lockout and the stroke indicator will display a lock symbol; in this state, the instrument will not allow the anvil release button to reopen the instrument jaws.
- The instrument can only achieve a maximum articulation angle of 45°. When using body structures or organs as a grounding surface, particular attention should be placed to the visual cues and tactile feedback received from the instrument. When the maximum angle is reached, the force will increase indicating the maximum angle has been

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reached. Avoid applying excessive pressure to the tissue as tissue damage or tissue trauma may occur.

- When dividing major vascular structures, be sure to adhere to the basic surgical principle of proximal and distal control.
- Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line.
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.
- Ensure tissue has not extended (extruded) proximal to the proximal black line on the instrument. Tissue forced into the instrument proximal to the black line may be transected without staples.
- If the closing trigger (1) is difficult to lock, reposition the instrument and take a smaller amount of tissue. Ensure that the proper reload selection has been made. (Refer to the ECHELON Reload Product Codes table.)
- Selection of the appropriate staple cartridge should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to fire and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.
- Attempting to force the trigger to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in increased force to fire or instrument failure.
- If the instrument locks out, the stroke count indicator will display a lockout symbol. Stop and push the red manual knife reverse switch downward to reverse the knife motion; the knife direction indicator will display an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. Squeeze firing trigger (2) completely until it rests on the closing trigger. The stroke count indicator will display “0” to indicate the knife has returned to its home position. Then press the anvil release button to release the jaws

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from the tissue. Close the instrument jaws by pressing the closing trigger (1), remove the instrument, and replace the cartridge. Firing through the lockout mechanism will break the instrument.

- If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the retracted position by verifying that the stroke count indicator displays “0” and knife direction indicator points towards the proximal side of the instrument, or that the knife blade indicator is in the home position. If the stroke count indicator or knife blade indicator is not in the home position, push the red manual knife reverse switch downward to reverse the knife motion and squeeze the firing trigger (2) completely until it rests on the closing trigger (1). Press the anvil release button. If the jaws do not open at this point, then gently pull the closing trigger (1) upward (away from the handle) until both firing and closing triggers return to their original positions
- If the firing mechanism becomes inoperative, do not continue to use the instrument.
- If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use the instrument.
- Before removing the instrument, be sure tissue is cleared from the jaws and then close the jaws.
- After removing the instrument, examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with electrocautery, manual sutures or other appropriate techniques.
- The instrument may be reloaded during a single procedure. Do not load the instrument more than 12 times for a total of 12 firings per instrument.
- Once the knife blade has been returned to its start position either manually or automatically, the cartridge is spent and must be discarded. If firing is attempted with a spent cartridge, the instrument will lock out. If resistance is felt, lockout has occurred. If the instrument locks out, the stroke count indicator will display a lockout symbol. Stop and push the red manual knife reverse switch downward to reverse the knife motion; the knife direction indicator will display an arrow pointing towards the proximal end of the instrument to indicate

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the knife is in the return mode. Squeeze firing trigger (2) completely until it rests on the closing trigger (1). The stroke count indicator will display "0" to indicate the knife has returned to its home position. Then press the anvil release button to release the jaws from the tissue. Close the instrument jaws by pressing the closing trigger (1), remove the instrument, and replace the cartridge. Firing through the lockout mechanism will break the instrument.

- When selecting the cartridge, careful consideration should be given to existing pathologic conditions as well as any pre-surgical treatment, such as radiotherapy, that the patient may have undergone. Certain conditions or preoperative treatments may cause changes in tissue thickness that would exceed the indicated range of tissue thickness for the standard choice of cartridge.

- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.

- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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How Supplied

The ECHELON families of endoscopic linear cutters (articulating and straight) are supplied sterile for single patient use. Discard after use.

The reloads are supplied sterile for single patient use. Discard after use.

The contents of cardboard packaging (boxes):

- Reload – 12 items.
- Instruction manual – 1 item (is supplied together with group package).

Note: for an explanation of the symbols on the package see the end of this manual.

Disposal

Disposal after use and expiring date (sterility term) must be carried out in accordance with accepted standards of medical practice and applicable local, state and federal laws and regulations.

After use, the instrument is a potential biohazard because after use it can be a potential biological hazard due to contamination with blood and other biological materials.

Disposal must be carried out in accordance with accepted medical practice standards, and relevant local, state, and Federal legislative acts and regulations according to SanPiN 2.1.7.2790-10, class B.

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2018 year / г.

Sarah Ward
23 Feb 2018
SARAH WARD
NOTARY PUBLIC
STATE OF OHIO
Comm. Expires
12-02-2022

ETHICON ECHELON FLEX

Reload ECHELON 45 mm, 6 rows:
black, white, blue, gold, green

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Please read all information carefully.


Failure to properly follow the instructions may lead to serious surgical consequences, such as leakage or disruption.

Important: This package insert is designed to provide instructions for use of the ECHELON ENDOPATH 45 mm Endoscopic Linear Cutter Reload with Gripping Surface Technology. It is not a reference to surgical techniques.

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
Instruction for use

Designation

ECHELON ENDOPATH 45 mm Endoscopic Linear Cutter (further Instrument) and Reload with Gripping Surface Technology (further reloads) are intended for transection, resection, and/or creation of anastomoses.

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Indications

The Instrument and reloads are intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures. They can be used with staple line or tissue buttressing materials. The instruments may also be used for transection and resection of liver parenchyma (hepatic vasculature and biliary structures), pancreas, kidney and spleen.

Area of application

It is used only in medical and preventive health care institutions.

Users

The instrument is intended for use only by highly qualified physicians with sufficient knowledge in the field of surgical techniques.

Contraindications

- Do not use the instruments on ischemic or necrotic tissue.
- Do not use the instruments on the aorta.
- Do not use any linear cutter on major vessels without making provision for proximal and distal control.
- Tissue thickness should be carefully evaluated before firing any stapler. Refer to the Reload Product Codes Table below for tissue compression requirement (closed staple height) for each staple size. If tissue cannot comfortably compress to the closed staple height listed in the table, or easily compresses to less than the closed staple height listed in the table, the tissue is contraindicated as it may be too thick or too thin for the selected staple size.
- These instruments are not intended for use when surgical stapling is contraindicated.

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Undesirable Side Effects / Residual Risks

Undesirable side effects and risks associated with powered surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

Description

Instrument and reloads are sterile, single patient use instruments that simultaneously cut and staple tissue. There are six staggered rows of staples, three on either side of the cut line. The 45 mm instruments have a staple line that is approximately 45 mm long and a cut line that is approximately 42 mm long. The shaft can rotate freely in both directions and an articulation mechanism on articulating instruments enables bending the distal portion of the shaft to facilitate lateral access of the operative site.

The instruments are shipped without a reload and must be loaded prior to use. A staple retaining cap on the reload protects the staple leg points during shipping and transportation. The instruments' lock-out feature is designed to prevent a used reload from being refired.

CAUTION: Do not load reloads into instruments more than 12 times for a maximum of 12 firings per instrument. The use of the instrument with staple line reinforcement material may reduce the number of firings.

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Reload Product Codes Chart

WARNING: Reloads may only be used with the ECHELON family of 45 mm Endoscopic Linear Cutters.

Reload Product Code	Open Staple Height	Closed Staple Height	Staple Line Length	Reload Color	Number of Staples	Staple Rows
GST45W	2.6 mm	1.0 mm	45 mm	White	70	6
GST45B	3.6 mm	1.5 mm	45 mm	Blue	70	6
GST45D	3.8 mm	1.8 mm	45 mm	Gold	70	6
GST45G	4.1 mm	2.0 mm	45 mm	Green	70	6
GST45T	4.2 mm	2.3 mm	45 mm	Black	70	6

MR Conditional

Non-clinical testing has demonstrated the staples made of titanium (Ti3Al2.5V) alloy are MR Conditional. A patient with the implanted staples can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 8,900 gauss/cm (89 T/m.)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)

Under the scan conditions defined above, the staples are expected to produce a maximum temperature rise of less than 2 °C after 15 minutes of continuous scanning.

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In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the staple line when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

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Operation Instructions

Verify compatibility of all instruments and accessories prior to using the instrument (refer to Warnings and Precautions).

Loading the Instrument

1 Using sterile technique, remove the instrument and reload from their respective packages. To avoid damage, do not flip the instrument or reload into the sterile field.

2 Prior to loading, ensure the instrument is in the open position (Illustration 2).

3 Examine the reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

CAUTION: Selection of the appropriate staple reload should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material.

WARNING: Loading the incorrect reload for the instrument size or model may result in transecting the tissue without stapling (e.g. loading a 45 mm reload into a 60 mm instrument).

4 Insert the new reload by sliding it against the bottom of the reload jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now loaded and ready for use (Illustration 3).

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The staple retaining cap ensures proper staple orientation and protects the staple leg points during shipping and transportation.

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5 After the device has been loaded, refer to the endoscopic linear cutters inserts for instructions for use.

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3 Push upward (toward the anvil) to unsnap the reload from the reload jaw. Discard the used reload (Illustration 4).

WARNING: Prior to reloading the instrument, rinse the anvil and reload jaw in sterile solution and then wipe the anvil and reload jaw to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and reload jaw.

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5 Insert the new reload by sliding it against the bottom of the reload jaw until the reload alignment tab

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stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use (Illustration 3).

CAUTION: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)

6 After the device has been loaded, refer to the endoscopic linear cutters inserts for instructions for use.

Standard Conventions Used

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a Warning or Caution statement. These statements are found throughout the documentation. These statements should be read before continuing to the next step in a procedure.

WARNING: A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.

CAUTION: A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.

Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.

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- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.

- When using other technologies (e.g., electrocautery) in the procedure, observe the precautions suggested by the original equipment manufacturer to avoid the hazards associated with their use.

- When selecting an instrument and reload, careful consideration should be given to existing pathologic conditions as well as any pre surgical treatment, such as radiotherapy, that the patient may have undergone. Certain conditions or preoperative treatments may cause changes in tissue thickness that would exceed the indicated range of tissue thickness for the standard choice of reload.

- Selection of the appropriate staple reload should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to fire and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.

- Inspect the package for shipping damage. Do not use an instrument, battery pack, or reload that has shipping damage.

- The instrument may be reloaded during a single procedure. Do not load the instrument more than 12 times for a maximum of 12 firings per instrument.

- ECHELON ENDOPATH™ 45 mm reloads may only be used with the ECHELON and ECHELON FLEX™ families of 45 mm endoscopic linear cutters.

- Loading the incorrect reload for the instrument size or model may result in transecting the tissue without

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stapling (e.g. loading a 45 mm reload into a 60 mm instrument).

- Prior to reloading the instrument, rinse the anvil and reload jaw in sterile solution and then wipe the anvil and reload jaw to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and reload jaw.

- Before using, remove the staple retaining cap and observe the surface of the reload. The reload must be replaced with another reload if any colored drivers are visible (if colored drivers are visible, the reload may not contain staples).

- When dividing major vascular structures, be sure to adhere to the basic surgical principle of proximal and distal control.

- Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line.

- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.

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- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single use devices may create a risk of contamination and/or cause infection or cross-infection, including but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Storage Conditions

Temperature: -20 °C – +25 °C
Relative Humidity: 10% – 65%

Transport Conditions

Temperature: -22 °C – +65 °C
Relative Humidity: 10% – 65%

Operating Environment Conditions

Temperature: +10 °C – +40 °C
Relative Humidity: 30% – 75%

Disposal

After expiration date of instrument and reloads (period of preservation of sterility) and/or in damaged packaging it must be disposed of as medical waste of class A according

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to SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" in accordance with the rules of the medical institution. After use, the product is a potential biological hazard from contamination with blood and other biological materials. Recycling should also be in accordance with accepted medical practice standards, and relevant local, state, and Federal legislative acts and instructions, as the waste classes B according to SanPiN 2.1.7.2790-10.

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Sterilization

Sterilization is made by gamma irradiation.

Sterility term – 3 years

Shelf life

3 years

Shelf for implantation: not limited

Standard configuration

The reloads is supplied sterile for single patient use. Discard after use.

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- reloads, 12 pieces.

Instruction manual (1 piece) is supplied together with group package.

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475 Calle C, Guaynabo, PR 00969 USA

Manufacturer warranty

A storage warranty period - the manufacturer guarantees the preservation of the sterility of the product and its suitability for the intended use, subject to the right

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transportation and storage within 3 years from the date of production.

Manufacturer is not liable for any damage or accidents resulting from improper use, human error or external emergencies. Regarding the product quality, please contact an authorized manufacturer's representative.

Authorized manufacturer's representative

LLC Johnson & Johnson
Krylatskaya 17, Building 2, 121614 Moscow, Russia
tel. +7(495)5807777, fax +7(495)5807878

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	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Manufacturer
	Catalogue Number
	Lot
	Authorized Representative in the USA
	Authorized Representative in the European Community
	Use Until Date
	Product passed procedure of declaration of conformity

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

Important: This package insert is designed to provide instructions for use of the ENDOPATH ECHELON 60 mm Endoscopic Linear Cutter Reloads with Gripping Surface Technology. It is not a reference to surgical techniques.

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Instruction for use

Designation

ECHELON ENDOPATH 60 mm Endoscopic Linear Cutter (further Instrument) and Reload with Gripping Surface Technology (further reloads) are intended for transection, resection, and/or creation of anastomoses.

Indications

The Instrument and reloads are intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures. They can be used with staple line or tissue buttressing materials. The instruments may also be used for transection and resection of liver parenchyma (hepatic vasculature and biliary structures), pancreas, kidney and spleen

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Contraindications

- Do not use the instruments on ischemic or necrotic tissue.
- Do not use the instruments on the aorta.
- Do not use any linear cutter on major vessels without making provision for proximal and distal control.

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- Tissue thickness should be carefully evaluated before firing any instrument. Refer to the Reload Product Codes Table below for tissue compression requirement (closed staple height) for each staple size. If tissue cannot comfortably compress to the closed staple height listed in the table, or easily compresses to less than the closed staple height listed in the table, the tissue is contraindicated as it may be too thick or too thin for the selected staple size.
- These instruments are not intended for use when surgical stapling is contraindicated.

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Device Description

The Instrument and reloads are sterile, single patient use instruments that simultaneously cut and staple tissue. There are six staggered rows of staples, three on either side of the cut line. The instruments have a staple line that is approximately 60 mm long and a cut line that is approximately 57 mm long. The shaft can rotate freely in both directions and an articulation mechanism on

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articulating instruments enables bending the distal portion of the shaft to facilitate lateral access of the operative site.

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Caution: Do not load instruments more than 12 times for a maximum of 12 firings per instrument. The use of the instrument with staple-line reinforcement material may reduce the number of firings.

Reload Product Codes Chart

WARNING: Reloads may only be used with the ECHELON family of 60 mm Endoscopic Linear Cutters.

Reload	Open Staple Height	Closed Staple Height	Staple Line Length	Reload Color	Number of Staples	Staple Rows
GST60W	2.6 mm	1.0 mm	60 mm	White	88	6
GST60B	3.6 mm	1.5 mm	60 mm	Blue	88	6
GST60D	3.8 mm	1.8 mm	60 mm	Gold	88	6
GST60G	4.1 mm	2.0 mm	60 mm	Green	88	6
GST60T	4.2 mm	2.3 mm	60 mm	Black	88	6

MR Conditional

Non-clinical testing has demonstrated the implantable staple made of titanium (Ti3Al2.5V) alloy in the ENDOPATH ECHELON 60 mm Endoscopic Linear Cutter Reload with Gripping Surface Technology reloads

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is MR Conditional. A patient with the implanted staples can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less.
- Highest spatial magnetic gradient field of 7.5 Tesla/m.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.99 W/kg for 15 minutes of scanning (per pulse sequence).

MRI Related Heating

In non-clinical testing, staple lines produced a temperature rise of less than 2 °C using the following conditions:

At 3-Tesla (Signa Hdxt General Electric Medical Systems MR scanner, software version 15.0_M4_910.a), a maximum MR system-reported whole body averaged SAR of 2.0 W/kg.

15 minutes of continuous MR scanning (per pulse sequence) using transmit/receive RF body coil.

Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the ENDOPATH ECHELON 60 mm Endoscopic Linear Cutter Reload with Gripping Surface Technology staple line. Therefore, optimization of MR imaging parameters to compensate for the presence of these staples may be necessary.

The worst case signal void size for the longest staple line (60 mm) was:

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1. Staple Retaining Cap
2. Anvil Jaw
3. Channel Jaw
4. Reload Alignment Slot
5. Reload Alignment Tab

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Loading the Instrument

- 1 Using sterile technique, remove the instrument and reload from their respective packages. To avoid damage, do not flip the instrument or reload into the sterile field.
- 2 Prior to loading, ensure the instrument is in the open position (Illustration 2).
- 3 Examine the reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

Caution: Selection of the appropriate staple cartridge should be based upon the combined thicknesses of both the tissue and the staple-line reinforcement material.

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4 Insert the new reload by sliding it against the bottom of the cartridge jaw until the cartridge alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. (Note: The staple retaining cap ensures proper staple orientation and protects the staple leg points during shipping and transportation.) The instrument is now loaded and ready for use (Illustration 3).

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5 Refer to the Instrument inserts for instructions for using the instrument after it is loaded.

Reloading the Instrument

1 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

2 Prior to reloading, ensure the instrument is in the open position (Illustration 2).

3 Push upward (toward the anvil) to unsnap the reload from the cartridge jaw. Discard the used reload (Illustration 4).

WARNING: Prior to reloading the instrument, rinse the anvil and cartridge jaw in sterile

solution and then wipe the anvil and cartridge jaw to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and cartridge jaw.

4 Examine the new reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

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5 Insert the new reload by sliding it against the bottom of the cartridge jaw until the cartridge alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use (Illustration 3).

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6 Refer to the ECHELON FLEX families of endoscopic linear cutters inserts for instructions for using the instrument after it is loaded.

Standard Conventions Used

The Use of Caution, Warning, and Note Statements

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a Caution, a Warning, or a Note statement. These statements are found throughout the documentation.

These statements should be read before continuing to the next step in a procedure.

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Caution: A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.

Note: A Note statement indicates an operating or maintenance problem, practice, or condition that is necessary to accomplish a task efficiently.

Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- When using other technologies (e.g., electrocautery) in the procedure, observe the precautions suggested by the original equipment manufacturer to avoid the hazards associated with their use.
- Inspect the package for shipping damage. Do not use an instrument, battery pack, or reload that has shipping damage.
- The instrument may be reloaded during a single procedure. Do not load the instrument more than 12 times for a maximum of 12 firings per instrument.
- ECHELON 60 mm Reloads may only be used with the ECHELON family of 60 mm Endoscopic Linear Cutters.
- Loading the incorrect reload for the instrument size or model may result in transecting the tissue without

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sealing (e.g. loading a 60 mm reload into a 45 mm instrument).

- After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)

- Prior to reloading the instrument, rinse the anvil and cartridge jaw in sterile solution and then wipe the anvil and cartridge jaw to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and cartridge jaw.

- (Powered Articulating Instruments) If the instrument locks out, the motor will stop. Release the firing trigger and slide the knife reverse switch forward to return the knife to the home position. In this position, the instrument should be removed, opened, and reloaded in order to continue.

- (Articulating Instruments) Do not attempt to articulate by pressing the front of the jaws against the grounding surface as tissue damage or tissue trauma may occur.

- Before using, remove the staple retaining cap and observe the surface of the reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)

- (Straight instruments) When placing the instrument through the trocar or incision, avoid advancing the firing trigger (2) accidentally as this action might activate the lockout. If this occurs, the instrument will not allow the anvil release button to reopen the instrument jaws. To recover from this condition, remove the instrument, pull the manual release lever; then press the anvil release button and replace the cartridge.

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- (Powered Articulating Instruments) When placing the instrument through the trocar or incision, avoid inadvertently pulling the red firing trigger lock and the firing trigger. The instrument may be partially or completely fired and will need to be reloaded before using on tissue.
- (Manual Articulating Instruments) When placing the instrument through the trocar or incision, avoid advancing the firing trigger (2) accidentally. If this occurs, the instrument will lockout and the stroke indicator will display a lock symbol; in this state, the instrument will not allow the anvil release button to reopen the instrument jaws.
- (Articulating Instruments) The instrument can only achieve a maximum articulation angle of 45°. When using body structures or organs as a grounding surface, particular attention should be placed to the visual cues and tactile feedback received from the instrument. When the maximum angle is reached, the force will increase indicating the maximum angle has been reached. Avoid applying excessive pressure to the tissue as tissue damage or tissue trauma may occur.
- When dividing major vascular structures, be sure to adhere to the basic surgical principle of proximal and distal control.
- Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line.
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.
- Ensure tissue has not extended (extruded) proximal to the proximal black line on the instrument.

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Tissue forced into the instrument beyond the black line may be transected without staples.

- If the closing trigger is difficult to lock, reposition the instrument and take a smaller amount of tissue. Ensure that the proper reload selection has been made. (Refer to the Reload Product Codes Chart.)

- (Manual Straight Instruments) The Endoscopic Linear Cutter is a multiple stroke firing device. Three complete strokes of the firing trigger will advance the staple line a full 60 mm and automatically return the knife blade. The manual firing release lever will return the knife blade and end the staple line before the full 60 mm staple line is applied.

- Selection of the appropriate staple cartridge should be based upon the combined thicknesses of both the tissue and the staple-line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to fire and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.

- (Powered Articulating Instruments) Attempting to force the device to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in motor stall and the knife will stop. If this happens, release the firing trigger, slide the knife reverse switch forward, and remove and reload the instrument. Then, position the instrument around a smaller tissue section or use a more appropriate reload (refer to the Reload Product Codes Chart).

- (Manual Instruments) Attempting to force the trigger to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in increased force to fire or instrument failure.

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- (Manual Articulating Instruments) If the instrument locks out, the stroke count indicator will display a lockout symbol. Stop and push the red manual knife return button downward to reverse the knife motion; the knife direction indicator will display an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. Squeeze firing trigger (2) completely until it rests on the closing trigger. The stroke count indicator will display "0" to indicate the knife has returned to its home position. Then press the anvil release button to release the jaws from the tissue. Close the instrument jaws by pressing the closing trigger (1), remove the instrument, and replace the cartridge. Firing through the lockout mechanism will break the instrument.

- (Powered Articulating Instruments) If the jaws do not automatically open after the anvil release switch is pressed, first ensure that the knife is in the home position. The position of the knife can be determined by observing the knife blade indicator under the cartridge jaw. If the knife blade indicator is not in the home position or the position of the knife cannot be determined, slide the knife return switch to activate the motor and return the knife to home position. Try opening the jaws again using the anvil release switch. If the jaws do not open at this point, then gently pull the closing trigger (1) upward (away from the handle) until both firing and closing triggers return to their original positions.

- (Articulating Instruments) If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the retracted position by verifying that the stroke count indicator displays "0" and knife direction indicator points towards the proximal side of the instrument, or that the knife blade indicator is in the home position. If the stroke count indicator or knife blade indicator is not in the home position, push the red manual knife return button downward to reverse the knife motion and squeeze the firing trigger (2) completely

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until it rests on the closing trigger (1). Press the anvil release button. If the jaws do not open at this point, then gently pull the closing trigger (1) upward (away from the handle) until both firing and closing triggers return to their original positions.

- Crossing of staple lines may shorten the life of the instrument.

- (Straight Instruments) If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the retracted position by verifying that the 3-stroke indicator or knife blade indicator is in the starting position. If the 3-stroke indicator or knife blade indicator is not in the starting position, pull the manual firing release lever until the 3-stroke indicator or knife blade indicator returns to the starting position. If the jaws do not open at this point, then gently pull the closing trigger (1) upward (away from handle) until both firing and closing triggers return to their original positions.

- The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the instrument.

- If the firing mechanism becomes inoperative, do not continue to use the instrument.

- If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use the instrument.

- Before removing the instrument, be sure tissue is cleared from the jaws and then close the jaws.

- Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.

- (Articulating Instruments) Once the knife blade has been returned to its start position either manually or

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automatically, the cartridge is spent and must be discarded. If firing is attempted with a spent cartridge, the instrument will lock out. If resistance is felt, lockout has occurred. If the instrument locks out, the stroke count indicator will display a lockout symbol. Stop and push the red manual knife return button downward to reverse the knife motion; the knife direction indicator will display

an arrow pointing towards the proximal end of the instrument to indicate the knife is in the return mode. Squeeze firing trigger (2) completely until it rests on the closing trigger (1). The stroke count indicator will display "0" to indicate the knife has returned to its home position. Then press the anvil release button to release the jaws from the tissue. Close the instrument jaws by pressing the closing trigger (1), remove the instrument, and replace the cartridge. Firing through the lockout mechanism will break the instrument.

- (Straight instruments) Once the knife blade has been returned to its start position either manually or automatically, the cartridge is spent and must be discarded. If firing is attempted with a spent cartridge, the instrument will lock out. If resistance is felt, lockout has occurred. If lockout occurs, stop, remove the instrument, and pull the manual release lever; then press the anvil release button and replace the cartridge. Firing through the lockout mechanism will break the instrument.

- When selecting the cartridge, careful consideration should be given to existing pathologic conditions as well as any preoperative treatments may cause changes in tissue thickness that would exceed the indicated range of tissue thickness for the standard choice of cartridge.

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- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only. Multiple patient use may compromise the device integrity or create a risk of contamination that, in turn, may result in patient injury or illness.

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Relative Humidity: 30% – 75%

Disposal

After expiration date of instrument and reloads (period of preservation of sterility) and/or in damaged packaging it must be disposed of as medical waste of class A according to SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" in accordance with the rules of the medical institution.

After use, the product is a potential biological hazard from contamination with blood and other biological materials.

Recycling should also be in accordance with accepted medical practice standards, and relevant local, state, and Federal legislative acts and instructions, as the waste classes B according to SanPiN 2.1.7.2790-10.

Maintenance and repair

This item is not applicable for the medical device (instrument and reloads).

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Sterilization

Sterilization is made by gamma irradiation.
Sterility term – 3 years

Shelf life

3 years
Shelf for implantation: not limited

Standard configuration

The reloads is supplied sterile for single patient use.
Discard after use.

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475 Calle C, Guaynabo, PR 00969 USA

Manufacturer warranty

A storage warranty period - the manufacturer guarantees the preservation of the sterility of the product and its suitability for the intended use, subject to the right transportation and storage within 3 years from the date of production.

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LLC Johnson & Johnson
Krylatskaya 17, Building 2, 121614 Moscow, Russia
tel. (495)5807777, fax (495)5807878

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Interpretation of symbols

	See Instructions For Use
	Single Patient Use
	Sterilized by Irradiation. Sterility Guaranteed Unless Package Opened or Damaged. Do Not Resterilize.
	Temperature limitation
	Humidity limitation
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Manufacturer/Date of Manufacture
	Catalogue Number
	Lot
	Authorized Representative in the USA
	Authorized Representative in the European Community

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	Use Until Date
	Compliance with requirements of European Directive 93/42/EES
	Product passed procedure of declaration of conformity
	Keep from moisture, store in a dry place
	Avoid of exposure direct sunlight. Keep away from sources of heat
	MR Conditional Compatibilité avec l'IRM (Imagerie par Résonance Magnétique) sous certaines conditions
	Manufacturer financially contributes to the cost of recovery and recycling

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М.П. / Stamp

INSTRUCTION FOR USE OF THE MEDICAL DEVICE

Selectable Cartridge/reloads for Linear Cutter

manufactured by
Ethicon Endo-Surgery, LLC,
475 Calle C, Guaynabo, Puerto Rico 00969, USA

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
2018 year / г.


Sarah Ward
16 March 2018
SARAH WARD
NOTARY PUBLIC
STATE OF OHIO
Comm. Expires
12-02-2022

Selectable Cartridge/reloads for Linear Cutter

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Important: This package insert is designed to provide instructions for use of the Straight Endoscopic Suturing Apparatus ECHELON 45 with knife and the staple cartridges. It is not a reference to surgical techniques.









ECHELON and Endopath are trademarks of Ethicon EndoSurgery.

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





Device Description

The Linear Cutter 55/75 mm is a sterile, single patient use instrument used in open surgical procedures. The instrument delivers six staggered rows of staples, three on either side of the cut line. The staple height is selectable to accommodate various tissue thicknesses using one cartridge/reload. The selectable staple height feature provides the ability to use one cartridge/reload to form staples with closed staple heights of 1.5 mm, 1.8 mm, or 2.0 mm.

The instrument has a safety lock-out feature that is designed to prevent it from firing if 1) no cartridge/ reload is loaded or 2) a used cartridge/reload is loaded. The instrument is shipped without a cartridge/reload loaded into it. The cartridge/reload must be loaded prior to use.

The cartridge/reload contains six rows of staples with the knife integrated in the cartridge/reload.

A staple retaining cap on the cartridge/reload protects the staples during shipping and transportation.

The 55 mm instrument creates a 61 mm staple line and a 58 mm cut line.

The 75 mm instrument creates a 81 mm staple line and a 78 mm cut line.

MR Conditional

Non-clinical testing has demonstrated the implantable staple made of titanium (Ti3Al2.5V) alloy in these devices are MR Conditional. A patient can be scanned safely immediately after placement of these staples, under the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Spatial gradient field of 720 Gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 2.7 W/kg for 15 minutes of scanning.

MRI Related Heating

In non-clinical testing, the implantable staple in the Linear Cutters produced a temperature rise of less than 2°C, at a maximum whole body averaged specific absorption rate (SAR) of 2.7 W/kg, as

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assessed by calorimetry for 15 minutes of MR scanning in a (3-Tesla/128-MHz, Excite, Software HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR scanner.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Proximate® Linear Cutters implantable staple. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) per staple extends approximately 3-mm relative to the size and shape of the implantable staple.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular
Signal Void Size (mm ²)	17	5	42	24

Linear Cutter 55/75 mm and Selectable Cartridges/Reloads Product Codes Chart

The product codes for the instruments and staple cartridges/reloads are as follows:

Instrument	Cartridge/Reload	Description	Number of Staples	Closed Staple Height and Selector Setting	Cartridge/Reload Color
NTLC55	SR55	55 mm Length	88	1.5 mm (Blue) 1.8 mm (Gold) 2.0 mm (Green)	Black
NTLC75	SR75	75 mm Length	118	1.5 mm (Blue) 1.8 mm (Gold) 2.0 mm (Green)	Black

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The instrument may be reloaded during a single procedure. Do not load the instrument more than twelve times for a maximum of twelve firings per instrument. The use of the instrument with staple line buttressing materials may reduce the number of firings.

Note: The Linear Cutter 55/75 mm is only compatible with the Linear Cutter 55/75 mm selectable cartridges/reloads.

Illustration and Nomenclature

- | | |
|----------------------------|---------------------------|
| 1. Anvil Half | 8. Staple Height Selector |
| 2. Locking Rib | 9. Cartridge/Reload Half |
| 3. Anvil Fork | 10. Firing Knob |
| 4. Staple Retaining Cap | 11. Gripping Surface |
| 5. Alignment/Locking Lever | 12. Alignment Tab |
| 6. Alignment Slot | 13. Gripping Surface |
| 7. Cartridge/Reload Fork | |

Instructions for Use

Verify compatibility of all instruments and accessories prior to use (refer to **Warnings and Precautions**).

- Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the cartridge/reload into the sterile field.
- If instrument is not separated into halves, separate instrument by completely disengaging the alignment/locking lever. (Illustration 2)
- Load the instrument by inserting the selectable cartridge/reload by placing the alignment tabs into the alignment slots and pivoting the selectable cartridge/reload onto the cartridge/reload fork. Snap the cartridge/reload into position. (Illustration 3)
- Remove staple retaining cap by grasping the edge of the staple retaining cap and lift straight up from the cartridge/reload. Discard the staple retaining cap. (Illustration 3)

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Note: The staple retaining cap ensures proper staple orientation and protects the staple leg points during shipping and transportation.

- 5 Prior to placing instrument across the tissue, select the desired closed staple height. In order to select the desired staple height, push the staple height selector so that it is aligned with desired color indicator, Blue 1.5 mm, Gold 1.8 mm or Green 2.0 mm. (Illustration 13) If tissue cannot be comfortably compressed to the closed staple height, or easily compresses to less than the closed staple height, either the setting or the tissue is contraindicated for the selected staple size or setting.

Caution: Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may require modification of the normally selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.

Caution: Tissue thickness should be carefully evaluated before firing any stapler.

- 6 Place the instrument across the tissue for transection or into the lumen to form an anastomosis. (Illustration 4)
 - 7 With the alignment/locking lever in the completely opened position, join the instrument halves together by aligning from either the front, center, or back of the instrument. (Illustration 5)
 - 8 To adjust tissue on the forks before firing, move the alignment/locking lever to the intermediate position. This allows maneuvering of the tissue while the instrument halves are joined.
- Note:** This step is optional. The alignment/locking lever may be moved from the completely opened to the completely closed position.

Caution: Ensure that the tissue lies flat between the forks. Any “bunching” of tissue along the reload or scale may result in an incomplete staple line. Tissue to be transected must be located between the arrows marked on the instrument jaw (Illustration 6). Any tissue located outside of the arrows is out of the stapling range.

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Caution: When positioning the device on the application site, ensure that no obstructions such as clips, stents, guide wires, and etc., are within the instrument anvils. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

Caution: Before firing, ensure that the cartridge/reload fork and the anvil fork are aligned.

- 9 Close the alignment/locking lever completely when the tissue is properly in place. (Illustration 7)

Caution: Attempting to force the locking lever to complete the closing stroke with too much tissue or thickened tissue may result in poor staple formation with the loss of staple line integrity and subsequent leakage, disruption, or poor healing. In addition, instrument damage or failure may result.

Caution: Unusually high closure force is an indication to open the instrument and inspect for tissue anomalies and hard objects or consider replacing the instrument.

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Note: The use of staple line buttressing materials with the instrument may require an increased force to close.

- 10 With the instrument closed, the firing knob is rotated to either side of instrument. (Illustration 12)

Note: In its pre-firing position, the firing knob can not be rotated from its pre-firing position unless the alignment/latching lever is engaged.

Caution: Do not forcibly move the firing knob when it is in the pre-firing position and/or the alignment/latching lever is not engaged, as this action may prevent the instrument from firing properly. The firing knob can not be rotated from its pre-firing position unless the alignment/latching lever is engaged.

Caution: When using a tissue or staple line buttressing material, the instructions of the manufacturer of the buttress material should be followed. The use of staple line buttressing materials with the instrument may require an increased force to fire and may reduce the number of firings.

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Caution: Selection of the appropriate staple cartridge/reload should be based upon the combined thicknesses of both the tissue and the staple line buttressing materials.

- 11 To fire the linear cutter, place the thumb on the firing knob and two fingers on the shoulders of the linear cutter. Fire the instrument by pushing the firing knob completely forward. (Illustration 8)

Note: The use of staple line buttressing materials with the instrument may require an increased force to fire and may reduce the number of firings.

Note: Crossing of staple lines may shorten the life of the instrument.

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- 12 Completely return the firing knob to the original pre-firing, "RETURN KNOB HERE" position. (Illustration 9)
- 13 Separate the instrument halves by opening the alignment/locking lever and remove the instrument from the patient. (Illustration 4)

Caution: After removing the instrument, examine the staple lines for hemostasis/pneumostasis and proper staple closure. Minor bleeding can be controlled with electrocautery, manual sutures or other appropriate techniques.

Loading and Reloading the Instrument

1. Using sterile technique, remove the cartridge/reload from the package. To avoid damage, do not flip the instrument into the sterile field.
2. If Firing Knob is not in "RETURN KNOB HERE" position, return the firing knob to the original prefiring, "RETURN KNOB HERE" position. (Illustration 9)

Note: The cartridge/reload cannot be inserted unless the firing knob is in its original position.

3. Separate the instrument halves by pulling open the alignment/locking lever.
4. Pull upward on the gripping surface and unsnap the used cartridge/reload from the cartridge/reload fork. Discard the used cartridge/reload. (Illustration 10)

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Caution: Prior to reloading the instrument, rinse the anvil and cartridge/reload jaw in sterile solution and then wipe the anvil fork and cartridge/reload fork to clear any formed but unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and cartridge/reload jaw.

5. Examine the new cartridge/reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the cartridge/reload.
6. Insert the new cartridge/reload by placing the alignment tabs into the alignment slots and pivoting the cartridge/reload into the cartridge/reload fork. Snap the cartridge/reload into position. Remove the staple retaining cap and discard. (Illustration 11)

Caution: After reloading, observe the surface of the new cartridge/reload. If the colored staple drivers are visible, replace with another reload. The staple drivers are used to push the staples out of the pocket, through the tissue to be stapled and into the forming anvil.

7. The instrument is now reloaded and ready for use.

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Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.
- Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may require modification of the normally selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.
- Tissue thickness should be carefully evaluated before firing any stapler.
- When dividing major vascular structures, be sure to adhere to the basic surgical principle of proximal and distal control.
- Remove the staple retaining cap before firing the instrument.
- Ensure that the knife is not exposed by returning the firing knob to the original "RETURN KNOB HERE" position.
- Before firing, ensure that the instrument cartridge/reload fork and anvil fork are aligned and the tissue is properly positioned.
- Ensure that the tissue lies flat between the forks. Any "bunching" of tissue along the reload or scale may result in an incomplete staple line. Tissue to be transected must be located between the

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- arrows marked on the instrument jaw (Illustration 6). Any tissue located outside of the arrows is out of the stapling range.
- The firing stroke must be completed. The instrument cannot be opened unless the firing knob is returned to its original “RETURN KNOB HERE” position.
- Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.
- Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple formation with the loss of staple line integrity and subsequent leakage, disruption, or poor healing. In addition, instrument damage or failure may result.
- Before closing the forks, be sure that tissue is cleared from the forks.
- After removing the instrument, examine the staple lines for hemostasis/pneumostasis and proper staple closure. Minor bleeding can be controlled with electrocautery, manual sutures or other appropriate techniques.
- Ensure that the firing knob is in the original “RETURN KNOB HERE” position during and after reloading.
- When positioning the device on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are within the instrument anvils. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.
- Crossing of staple lines may shorten the life of the instrument.
- When using a tissue or staple line buttressing material, the instructions of the manufacturer of the buttress material should be followed. The use of staple line buttressing materials with the instrument may require an increased force to fire and may reduce the number of firings.
- Selection of the appropriate staple height should be based upon the combined thicknesses of both the tissue and the staple line buttressing materials.
- When selecting the staple height, using the selector lever (Illustration 13), careful consideration should be given to existing

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Reloads for Curved Cutter Stapler CONTOUR

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2018 year / г.


Sarah Ward
10 March 2018
SARAH WARD
NOTARY PUBLIC
STATE OF OHIO
Comm. Expires
12-02-2022

Reloads for Curved Cutter Stapler CONTOUR

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contour[®]
Curved Cutter Stapler and Reloads

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


Important: This package insert is designed to provide instructions for use of the CONTOUR[™] Curved Cutter Stapler and Reloads. It is not a reference to surgical stapling techniques.

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Instructions for use

This IFU applies to

Reloads for Curved Cutter Stapler CONTOUR, versions:

Blue reload for tissue of thickness 1.5 mm **CR40B**

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Instrument		of Staples		Color	Thickness
CS40B	CR40B	Standard	46	Blue	1.5 mm
CS40G	CR40G	Thick	46	Green	2.0 mm

NOTE: The blue (standard) and green (thick) reloads can be used interchangeably with the instruments.

MR Conditional

Non-clinical testing has demonstrated the implantable staple made of titanium (Ti3Al2.5V) alloy in these devices is MR Conditional. A patient with the implanted staples can be scanned safely immediately after placement of these staples, under the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Spatial gradient field of 720 Gauss/cm.
- Maximum whole body averaged specific absorption rate (SAR) of 2.7 W/kg for 15 minutes of scanning.

MRI Related Heating

In non-clinical testing, the implantable staple made of titanium (Ti3Al2.5V) alloy in the Contour® Curved Cutter Staplers produced a temperature rise of less than 2°C at a maximum whole body averaged specific absorption rate (SAR) of 2.7 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a (3-Tesla/128-MHz, Excite, Software HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR scanner.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Contour® Curved Cutter Staplers implantable staple. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) per staple extends approximately 3-mm relative to the size and shape of the implantable staple.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

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Illustration and Nomenclature (Illustration 1)

- | | |
|---------------------------|----------------------|
| 1. Curved Head | 7. Handle |
| 2. Cartridge Module | 8. Closure Trigger |
| 3. Retaining Pin | 9. Firing Trigger |
| 4. Push Rod | 10. Cartridge Module |
| 5. Retaining Pin Actuator | 11. Cartridge Anvil |
| 6. Release Button | 12. Staple Retainer |

General Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument.

(Refer to **Warnings and Precautions.**)

- Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- Remove the staple retainer from the instrument by pushing on the pad labeled "Push to Remove". (Illustration 2). Hold the cartridge when removing the retainer. Check that the cartridge is fully loaded into the instrument (press down) after the retainer has been removed. Discard the staple retainer.
- Grasp the instrument handle in the palm of the hand with the fingers around the closure trigger.
Position tissue to be cut and stapled between the anvil and the cartridge.
NOTE: Any tissue covering the hole in the anvil will be pierced by the retaining pin. Reposition the tissue if necessary.
- Either manually advance the retaining pin using the actuator button on top of the handle to capture the tissue within the jaw opening (Illustration 3A), or automatically advance the retaining pin by squeezing the closure trigger (Illustration 3B). Note that there is an audible click halfway through the closure stroke indicating that the device is in the intermediate position. In the intermediate position, the pin is fully seated in the anvil capturing the tissue, and the jaws are partially open. Reposition tissue within the instrument if desired. The closure trigger can be released at this point, and the instrument will

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maintain its jaw opening to allow for final positioning of the tissue to be cut and stapled. The retaining pin is fully engaged.

Caution: Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line. NOTE: If using the automatic advancement of the retaining pin, ensure that the tissue is correctly placed to avoid accidental piercing of tissue. NOTE: Do not grasp the firing trigger before the instrument is to be fired.

NOTE: Recommended operation is to manually advance the Retaining Pin.

- 5** Squeeze the closure trigger and handle together until the closure trigger is latched. Listen for an audible click. When the closure trigger is latched, the firing trigger moves to the ready-to-fire position (Illustration 4). The cartridge has clamped onto the tissue which is ready to be cut and stapled. Ensure that the closing stroke is completed. All fingers should be completely removed from the closing trigger to ensure that the closing system holds.

NOTE: If the tissue needs to be repositioned within the instrument before stapling, open the jaws by squeezing the handle and closure trigger together slightly, push the release button, and slowly release grasp of the closure trigger. The closure trigger will return to the fully open position and the jaws will release the tissue but the retaining pin will remain extended. Tissue can now be repositioned. NOTE: The Contour Curved Cutter Stapler is designed with a feature which allows partial closing if a used reload, no reload, or an incorrectly inserted reload is in the instrument. In these cases, the device will close approximately one-third of the stroke and stop, but will return to the full open position when the closure trigger is released. The used or misloaded cartridge module should then be replaced with a new reload or, if no cartridge was present, a cartridge should be loaded in the instrument. After following the above steps, any device which does not close either partially or completely should not be used.

NOTE: After closing, pausing for 15 seconds prior to firing may result in better tissue compression and better staple form.

NOTE: If closure trigger (1) is difficult to lock, reposition the instrument and take a smaller amount of tissue. (See Contraindications

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- 6 With the handle and closure trigger grasped in the palm of the hand, place all fingers around the firing trigger. Before firing, check to ensure the retaining pin is seated in the anvil.

Caution: Always ensure complete closing prior to firing.

Caution: Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

Caution: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

Caution: Care should be taken not to fire this instrument when the resulting staple line will be in tension.

- 7 Completely fire the instrument by squeezing the firing trigger until it touches the closure trigger (plastic to plastic), signifying that the instrument has fired the staples and knife blade, and the tissue has been transected. (Illustration 5) The firing trigger automatically returns to the intermediate position as soon as the hand releases it, leaving the closure trigger in the fully closed position. Make sure that the release button is not pressed during firing as proper formation of staples and cut line may be compromised.

Caution: Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.

Caution: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

- 8 After firing the instrument, open the jaws of the instrument by squeezing the closure trigger and pushing the release button with the thumb. For controlled release, while holding the release button down, open your hand to release the closure trigger. (Illustration 6) The cartridge drivers are exposed as a spent cartridge indicator. If the instrument becomes locked onto tissue and will not open by pressing

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the release button, the device can be opened by depressing the release button and pulling the closure trigger simultaneously.

Caution: Examine the staple line for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

Caution: Following a division of the rectum, it is recommended that a leak test is performed prior to the creation of the anastomosis.

- 9 Retract the retaining pin manually. Ensure that the staple line and the cut line are complete. Remove the instrument. (Illustration 7)

Muscle Stapling Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument.

(Refer to Warnings and Precautions.)

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 Remove the staple retainer from the instrument by pushing on the pad labeled "Push to Remove". (Illustration 2). Hold the cartridge when removing the retainer. Check that the cartridge is fully loaded into the instrument (press down) after the retainer has been removed. Discard the staple retainer.
- 3 Grasp the instrument handle in the palm of the hand with the fingers around the closure trigger. Position tissue to be cut and stapled between the anvil and the cartridge.
NOTE: Any tissue covering the hole in the anvil will be pierced by the retaining pin. Reposition the tissue if necessary.
- 4 Recommended operation, in skeletal muscle, is to manually advance the Push Rod using the Retaining Pin Actuator. Squeeze the closure trigger (Illustration 3B). Note that there is an audible click halfway through the closure stroke indicating that the device is in the intermediate position. In the intermediate position, the pin is fully seated in the anvil capturing the tissue, and the jaws are partially open. Reposition tissue within the instrument if desired. The closure trigger can be released at this point, and the instrument will maintain its jaw opening to allow for final positioning of the tissue to be cut and stapled. The retaining pin is fully engaged.

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Caution: Make sure tissue to be stapled is properly positioned in the jaws before closing the jaws or stapling. Bunching, tension, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

NOTE: In thick skeletal muscle tissue, two approaches are recommended for device closure:

- A single slow, controlled and continuous squeeze of the Closure Trigger • Progressively squeezing and relaxing (pumping) the Closure Trigger NOTE: Do not grasp the firing trigger before the instrument is to be fired.

- 5 Squeeze the closure trigger and handle together until the closure trigger is latched. Listen for an audible click. When the closure trigger is latched, the firing trigger moves to the ready-to-fire position (Illustration 4). The cartridge is clamped onto the tissue, which is ready to be stapled and cut. Ensure that the closing stroke is completed. All fingers should be completely removed from the closing trigger to ensure that the closing system holds.

NOTE: If the closure trigger is difficult to lock, reposition the instrument and take a smaller amount of tissue.

NOTE: Do not use excessive force to close the instrument. Forcing the instrument to fully close, with excess tissue, may result in the closure trigger breaking.

NOTE: If the tissue needs to be repositioned within the instrument before stapling, open the jaws by squeezing the handle and closure trigger together slightly, push the release button, and slowly release grasp of the closure trigger. The closure trigger will return to the fully open position and the jaws will release the tissue but the retaining pin will remain extended. Tissue can now be repositioned. NOTE: The muscle stapling Contour Curved Cutter Stapler is designed with a feature that allows partial closing if a used cartridge reload, no cartridge reload, or an incorrectly inserted cartridge reload is in the instrument. In these cases, the device will close approximately one-third of the stroke and stop, but will return to the full open position when the closure trigger is released. The used or misloaded cartridge module should then be replaced with a new cartridge reload or, if no cartridge was present, a cartridge should be loaded in the instrument. After following the above steps, any device that does not close either partially or completely should not be used.

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NOTE: After closing, pausing for 15 seconds prior to firing may result in better tissue compression and better staple form.

- 6 With the handle and closure trigger grasped in the palm of the hand, place all fingers around the firing trigger. Before firing, check to ensure the retaining pin is seated in the anvil.

Caution: Always ensure complete jaw closure prior to firing.

Caution: Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

Caution: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

Completely fire the instrument by squeezing the firing trigger until it touches the closure trigger (plastic to plastic), signifying that the instrument has fired the staples and knife blade, and the tissue has been transected. (Illustration 5) The firing trigger automatically returns to the intermediate position as soon as the hand releases it, leaving the closure trigger in the fully closed position. Make **sure that the release button is not pressed during firing as proper formation of staples and cut line may be compromised.**

Caution: Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.

Caution: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

- 7 After firing the instrument, open the jaws of the instrument by squeezing the closure trigger and pushing the release button with the thumb. For controlled release, while holding the release button down, open your hand to release the closure trigger. (Illustration 6) The cartridge drivers are exposed as a spent cartridge indicator. If the instrument becomes locked onto tissue and will not open by pressing the release button, the device can be opened by depressing the release button and pulling the closure trigger simultaneously.

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Caution: Examine the staple line for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

- 8 Retract the retaining pin manually. Ensure that the staple line and the cut line are complete. Remove the instrument. (Illustration 7)

Reloading the Contour Curved Cutter Stapler

- 1 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.
- 2 Depress the release button to ensure the instrument is in the open position. Ensure that the retaining pin has been retracted.
- 3 Remove the used reload from the instrument. Grasp the top of the reload and lift upward, unsnapping the reload from the jaws. (Illustration 8) Properly discard the used reload. Thoroughly clean the instrument by vertically submerging completely the anvil and cartridge channel in a sterile solution and swish vigorously. After swishing, visually inspect the anvil and cartridge channel and remove all residual staples and/or foreign matter from instrument prior to reloading.

Caution: Failure to properly clean the instrument prior to reloading could damage the instrument; compromise cut line and staple formation of subsequent firings.

- 4 Examine the new reload for the presence of a staple retainer. If the retainer is not in place, discard the reload. Do not remove the staple retainer until the reload cartridge is inserted into the device.

NOTE: The staple retainer cannot be removed until the cartridge is correctly loaded into the device.

The loaded device cannot be used until the staple retainer is removed.

- 5 Insert the new reload into the metal housing and snap into position. (Illustration 9) The tracks on each side of the reload should be used as guides to align the reload within the jaws of the instrument. When the reload is properly aligned, push the reload into the instrument until it is fully seated. (Illustration 10) Remove the staple retainer. (Illustration 11) Check that the reload is held firmly within the jaws.

NOTE: If the cartridge is removed from the device, whether the cartridge is spent or not, and if the staple retainer has already been removed from the cartridge, the cartridge cannot be reloaded into the device.

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6 After reloading and after removing the staple retainer, observe the surface of the new reload. If drivers protrude out of the reload, replace with another reload.

7 The stapler is now reloaded and ready for use.

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- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together during a procedure, verify compatibility prior to initiation of procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. **Do not** immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.
- Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.
- Do not fire device without a loaded cartridge.
- Be sure to remove the staple retainer prior to inserting tissue and closing the device.
- If the staple retainer is not in place in a cartridge module, discard the cartridge module and replace with a new one. Do not remove the staple retainer from the cartridge until the reload cartridge has been inserted into the instrument.
- Before closing the instrument on tissue, observe the surface of the new reload where the staples will emerge. If colored drivers protrude out of the reload, replace with another reload.
- Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.
- Ensure that unintentional tissue is not cut.
- Always ensure complete closing prior to firing.

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- Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.
- Make sure that the release button is not pressed during firing as proper formation of staples and cut line may be compromised.
- Care should be taken not to fire this instrument when the resulting staple line will be in tension.
- Before firing, check to ensure the retaining pin is seated in the anvil. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.
- Do not fire the instrument unless the closure trigger is properly latched against the handle.
- The firing trigger must be pulled back completely against the closure trigger to properly fire the instrument. Do not partially fire the instrument.
- Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.
- The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.
- After firing and release and before removing the instrument, check that the cut is complete.
- Retract the retaining pin manually prior to removing the device.
- Examine the staple lines for hemostasis and complete closure. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- Following a division of the rectum, it is recommended that a leak test is performed prior to the creation of the anastomosis.

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