**BD**

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TriPath Imaging Inc.,  
780 Plantation Drive Burlington,  
North Carolina 27215, USA

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
М.П. / Stamp

|                                          |                                             |
|------------------------------------------|---------------------------------------------|
| State of                                 | <u>NORTH CAROLINA</u>                       |
| County of                                | <u>DURHAM</u>                               |
| Sworn to and subscribed before me on the | <u>23</u> day of <u>MARCH</u> , <u>2016</u> |
| Notary Public's Signature                | <u>Ann Lilley</u>                           |
| My Commission Expires                    | <u>12/29/2018</u>                           |

## OPERATIONAL DOCUMENTATION

(instruction for use)

for the medical device:

### Reagents for preparation and coloring samples for cytological tests

Manufactured by TriPath Imaging Inc.,  
780 Plantation Drive Burlington, North Carolina 27215, USA

## 1. Name of the medical device

Reagents for preparation and coloring samples for cytological tests

1. BD PrepStain Cytology Stain kit (gynecological material) in plastic vials:
  - Hematoxylin Stain 0.75 – 1 vial, 480 mL
  - EA/OG Combo Stain – 1 vial, 480 mL
2. Alcohol Blend Rinse – plastic container, 1,7 L
3. BD SurePath Preservative Fluid - plastic container, 3,6 L
4. TRIS Buffered Saline
5. Density Reagent in a vial - 480 mL.
6. BD CytoRich Red Preservative - plastic container, 3,6 L.
7. BD CytoRich Blue Preservative - plastic container, 3,6 L.
8. BD PrepStain Non-GYN Stain kit (non-gynecological material) in plastic vials:
  - Hematoxylin Stain 0,5 - 1 vial, 480 mL;
  - EA/OG Combo Stain - 1 vial, 480 mL.
9. BD SurePath Preservative Fluid in vials 10 mL (kit - 25 vials).
10. BD SurePath Preservative Fluid in vials 10 mL (kit - 500 vials).
11. BD SureDetect Counterstains:
  - Hematoxylin Stain in a vial 18 mL;
  - Blue Stain in a vial 18 mL.
12. BD ProEx C Immunocytochemical Test:
  - vial 15 mL – primary antibodies to protein supported MCM2 and to TOP2 $\alpha$ .
13. BD SureDetect Detection Reagents:
  - vial 15 mL - Peroxidase Block reagent;
  - vial 15 mL – Protein Block reagent;
  - vial 15 mL – Secondary Mouse Probe antibodies;
  - vial 15 mL – Reagent polymer based on peroxidase;
  - vial 18 mL – DAB Substrate Buffer (solvent for 3,3'- diaminobenzidine);
  - vial 1.0 mL - DAB Chromagen (dye, 3,3'- diaminobenzidine).
14. BD SureDetect SiHa Cell Control in a vial 8 mL for preparation of samples.
15. BD SureDetect Slide Preparation Buffer (10x).
16. BD ProEx C Immunohistochemical Test in a vial 7 mL.

## 2. Place of manufacture of the medical device

Legal manufacturer:

TriPath Imaging Inc., 780 Plantation Drive Burlington, North Carolina 27215, USA

Helping all people  
live healthy lives

#### Manufacturing sites:

1. TriPath Imaging Inc., 780 Plantation Drive, Burlington, North Carolina, 27215, USA.
2. TriPath Imaging Inc., 4025 Stirrup Creek Drive, Suite 400, Durham, North Carolina, 27703, USA.
3. TriPath Imaging Inc., 1022 Corporate Park Drive, Mebane, North Carolina 27302, USA.

### 3. Medical device classification

The European Directive on In Vitro Diagnostic Medical Devices: 98/79/EC

Class of potential risk of use of the medical device: I.

GMDN code: 12955.

### 4. Type of contact with the human body:

No contact with the body.

### 5. Intended use of the medical device

Reagents for preparation and coloring samples for cytological tests are used to carry out a specific immunostaining and detection during cytological tests.

1. BD PrepStain Cytology Stain kit is used for staining gynecological material and consists of Hematoxylin Stain 0.75 and EA/OG Combo Stain. Both stains are standard stains in cytology. The BD Hematoxylin Stain (hematoxylin) is applied, aspirated and reapplied to stain. BD Hematoxylin is the first stain used in the staining process. The BD EA/OG Combo Stain (EA/OG) is applied, aspirated, and reapplied to stain. EA/OG is the second stain used in the staining process.
2. Alcohol Blend Rinse is used for washing samples during staining process.
3. BD SurePath Preservative Fluid is used to collect and preserve gynecologic samples.
4. TRIS Buffered Saline is used for washing of cytological samples during staining process.
5. Density Reagent is used for washing of cytological samples from additional mixtures.
6. BD CytoRich Red Preservative is intended to preserve cells and small tissue fragments in suspension.
7. BD CytoRich Blue Preservative is intended for preservation of cytology cells in suspension and their subsequent preparation for cytological evaluation.
8. BD PrepStain Non-GYN Stain kit is used for staining non-gynecological material and consists of Hematoxylin Stain 0,5 and EA/OG Combo Stain.

Helping all people  
live healthy lives

9. BD SurePath Preservative Fluid in vials 10 mL (kit - 25 vials) is used to collect and preserve gynecologic samples.
10. BD SurePath Preservative Fluid in vials 10 mL (kit - 500 vials) is used to collect and preserve gynecologic samples.
11. BD SureDetect Counterstains is intended for the qualitative evaluation of aberrant S-Phase induction in cervical cytology specimens collected in SurePath Preservative Fluid. The test results provide adjunctive information for cervical cytology diagnosis. BD SureDetect Counterstains consists of Hematoxylin Stain and Blue Stain.
12. BD ProEx C Immunocytochemical Test is intended for the qualitative evaluation of aberrant S-Phase induction in cervical cytology specimens collected in SurePath Preservative Fluid. BD ProEx C Immunocytochemical Test consists of primary antibodies to protein supported MCM2 and to TOP2 $\alpha$ .
13. BD SureDetect Detection Reagents are intended for the qualitative evaluation of aberrant S-Phase induction in cervical cytology specimens collected in SurePath Preservative Fluid. BD SureDetect Detection Reagents are Peroxidase Block reagent, Protein Block reagent, Secondary Mouse Probe antibodies, Reagent polymer based on peroxidase, DAB Substrate Buffer (solvent for 3,3'-diaminobenzidine), DAB Chromagen (dye, 3,3'-diaminobenzidine).
14. BD SureDetect SiHa Cell Control may be used in immunoassays that require an epithelial-based cell line as an assay staining control. This product is to be used in slide-based assays to allow the user to verify the staining procedure.
15. BD SureDetect Slide Preparation Buffer is to be used on cytology and histology specimens, which have been mounted on glass for target retrieval prior to immunocytochemistry (ICC) or immunohistochemistry (IHC) procedures. Cytology specimens should be liquid-based cytology specimens prepared on glass slides. Histology specimens should be formalin-fixed paraffin-embedded tissue sections mounted on glass slides.
16. BD ProEx C Immunohistochemical Test contains an antibody reagent for use in a three-step manual or automated immunohistochemical staining procedure in formalin-fixed paraffin-embedded tissue biopsy specimens.

#### **Indications for administration:**

The device is not purposed for direct application on the patient.

Reagents for preparation and coloring samples for cytological tests are used only for IN VITRO (cytological) diagnostics of gynecological and non-gynecological diseases.

#### **Contraindications:**

No.

## 6. Sphere of application

Laboratory IN VITRO diagnostics.

## 7. Description and technical specification of the device

Reagents for preparation and coloring samples for cytological tests:

1. BD PrepStain Cytology Stain kit consists of Hematoxylin Stain 0.75 and EA/OG Combo Stain. Both stains are standard stains in cytology. High Density Polyethylene (HDPE) bottles and caps, proprietary formula for the Hematoxylin and EA/OG combo stain.

- Hematoxylin Stain 0.75 – 1 vial, 480 mL

### Composition of Hematoxylin Stain 0.75

|                   |                    |       |
|-------------------|--------------------|-------|
| CAS: 107-21-1     | ethylene glycol    | 10,5% |
| EINECS: 203-473-3 | Xn R22             |       |
|                   | Acute Tox. 4, H302 |       |

### Physical and chemical properties of Hematoxylin Stain 0.75

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Violet                                        |
| Odour:                                  | Characteristic                                |
| pH-value                                | 2,6                                           |
| Flash point                             | >100 °C                                       |
| Ignition temperature                    | 410 °C                                        |
| Self-igniting                           | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Density:                                | 1.02 g/cm <sup>3</sup>                        |
| Solubility in / Miscibility with Water: | Soluble                                       |

- EA/OG Combo Stain – 1 vial, 480 mL

### Composition of EA/OG Combo Stain

|                                   |                                                                                                                                                  |       |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 64-17-5<br>EINECS: 200-578-6 | ethanol<br>F R10<br>Flam. Liq. 2, H225                                                                                                           | 69,7% |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 25,0% |
| CAS: 64-19-7<br>EINECS: 200-580-7 | acetic acid<br>C R35                                                                                                                             | 2,0%  |

Helping all people  
live healthy lives

|  |                                                |  |
|--|------------------------------------------------|--|
|  | R10<br>Flam. Liq. 3, H226; Skin Corr. 1A, H314 |  |
|--|------------------------------------------------|--|


### Physical and chemical properties of EA/OG Combo Stain

| Property                                | Comments                                                                                    |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                             |                                                                                             |
| Form:                                   | Liquid                                                                                      |
| Colour:                                 | Dark green                                                                                  |
| Odour:                                  | Alcohol-like                                                                                |
| pH-value                                | 5,1 +/- 0,2                                                                                 |
| Boiling point/Boiling range             | 73 °C                                                                                       |
| Flash point                             | 12 °C                                                                                       |
| Ignition temperature                    | 425 °C                                                                                      |
| Self-igniting                           | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 3,5 Vol %                                                                                   |
| Upper:                                  | 44,0 Vol %                                                                                  |
| Density:                                | 0,826 g/cm <sup>3</sup>                                                                     |
| Solubility in / Miscibility with Water: | Not miscible or difficult to mix                                                            |
| Solvent content:                        |                                                                                             |
| Organic solvents:                       | 96,7 %                                                                                      |
| Water:                                  | 2,5 %                                                                                       |

2. Alcohol Blend Rinse – plastic container, 1,7 L. High Density Polyethylene (HDPE) bottle and 38/400 pp White r/m f217 cap, proprietary formula for the alcohol blend rinse.

Helping all people  
live healthy lives

### Composition of Alcohol Blend Rinse

|                                   |                                                                                                                                                  |       |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 67-63-0<br>EINECS: 200-661-7 | isopropanol<br>Xi R36; F R11<br>R67<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                                   | 50,0% |
| CAS: 64-17-5<br>EINECS: 200-578-6 | ethanol<br>F R10<br>Flam. Liq. 2, H225                                                                                                           | 48,0% |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 2,0%  |

### Physical and chemical properties of Alcohol Blend Rinse


| Property                                | Comments                                                                                    |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                             |                                                                                             |
| Form:                                   | Liquid                                                                                      |
| Colour:                                 | Clear                                                                                       |
| Odour:                                  | Characteristic                                                                              |
| Boiling point/Boiling range             | 78 °C                                                                                       |
| Flash point                             | 10 °C                                                                                       |
| Ignition temperature                    | 425 °C                                                                                      |
| Self-igniting                           | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 2,0 Vol %                                                                                   |
| Upper:                                  | 15,0 Vol %                                                                                  |
| Density                                 | 0.79 g/cm <sup>3</sup>                                                                      |
| Solubility in / Miscibility with Water: | Soluble                                                                                     |
| Solvent content:                        |                                                                                             |
| Organic solvents:                       | 100,0 %                                                                                     |

Helping all people  
live healthy lives

3. BD SurePath Preservative Fluid - plastic container, 3,6 L. High Density Polyethylene (HDPE) bottle and 38/400 pp White r/m f217 cap, proprietary formula for the preservative fluid.

#### Composition of BD SurePath Preservative Fluid

|                                   |                                                                                                                                                  |       |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 64-17-5<br>EINECS: 200-578-6 | ethanol<br>F R10<br>Flam. Liq. 2, H225                                                                                                           | 21,7% |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 1,2%  |
| CAS: 67-63-0<br>EINECS: 200-661-7 | isopropanol<br>Xi R36; F R11<br>R67<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                                   | 1,1%  |

#### Physical and chemical properties of BD SurePath Preservative Fluid


| Property                     | Comments     |
|------------------------------|--------------|
| Appearance:                  |              |
| Form:                        | Liquid       |
| Colour:                      | Clear        |
|                              | Colourless   |
| Smell:                       | Alcohol-like |
| Freezing point:              | 7,2          |
| Boiling point/Boiling range: | 86 °C        |

Helping all people  
live healthy lives

|                                         |                                                                                             |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Flash point                             | 35 °C                                                                                       |
| Ignition temperature:                   | 425,0 °C                                                                                    |
| Self-igniting                           | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 3,5 Vol %                                                                                   |
| Upper:                                  | 15,0 Vol %                                                                                  |
| Vapour pressure at 20 °C:               | 59,0 hPa                                                                                    |
| Solubility in / Miscibility with Water: | Soluble                                                                                     |
| Solvent content:                        |                                                                                             |
| Organic solvents:                       | 24,2 %                                                                                      |
| Water:                                  | 74,7 %                                                                                      |
| Solids content:                         | 0,9 %                                                                                       |

4. TRIS Buffered Saline is used for washing of cytological samples during staining process.

#### Composition of TRIS Buffered Saline


|                   |                                                          |        |
|-------------------|----------------------------------------------------------|--------|
| CAS: 77-86-1      | tris base                                                | 17,38% |
| EINECS: 201-064-4 | Xi R36/38                                                |        |
|                   | Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335 |        |

#### Physical and chemical properties of TRIS Buffered Saline

| Property    | Comments       |
|-------------|----------------|
| Appearance: |                |
| Form:       | Solid          |
| Color:      | White          |
| Odor:       | Characteristic |

|                                         |                                               |
|-----------------------------------------|-----------------------------------------------|
| Self-igniting                           | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

5. Density Reagent in a vial - 480 mL. Density Reagent is a liquid on the ground of sodium azide (0.1%), deionized water (>90%), sodium chloride (<1%), hydroxyethyl starch (<10%) and it is used to clear the smears from admixtures.

### Composition of Density Reagent

|                                      |                                                                                                               |      |
|--------------------------------------|---------------------------------------------------------------------------------------------------------------|------|
| CAS: 26628-22-8<br>EINECS: 247-852-1 | sodium azide<br>T+ R28; N R50/53<br>R32<br>Acute Tox. 2, H300; Aquatic Acute 1, H400; Aquatic Chronic 1, H410 | 0,1% |
|--------------------------------------|---------------------------------------------------------------------------------------------------------------|------|

### Physical and chemical properties of Density Reagent

| Property    | Comments       |
|-------------|----------------|
| Appearance: |                |
| Form:       | Liquid         |
| Colour:     | Clear          |
| Odour:      | Characteristic |
| pH-value:   | 7,0            |
| Flash point | >100 °C        |


Helping all people  
live healthy lives

|                                         |                                               |
|-----------------------------------------|-----------------------------------------------|
| Self-igniting                           | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Density:                                | 1,024 g/cm <sup>3</sup>                       |
| Solubility in / Miscibility with Water: | Soluble                                       |

6. BD CytoRich Red Preservative - plastic container, 3,6 L. High Density Polyethylene (HDPE) bottle and 38/400 pp White r/m f217 cap, proprietary formula for the preservative fluid CytoRich red. BD CytoRich Red Preservative is an alcohol-based hemolytic fixative. It is for cytology and small biopsy processing. Its ingredients include alcohols, formaldehyde, and buffers.

#### Composition of BD CytoRich Red Preservative

|                                    |                                                                                                                                                  |       |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 67-63-0<br>EINECS: 200-661-7  | isopropanol<br>Xi R36; F R11<br>R67<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                                   | 23,3% |
| CAS: 67-56-1<br>EINECS: 200-659-6  | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 10,0% |
| CAS: 107-21-1<br>EINECS: 203-473-3 | ethylene glycol<br>Xn R22<br>Acute Tox. 4, H302                                                                                                  | 6,7%  |

#### Physical and chemical properties of BD CytoRich Red Preservative

| Property | Comments |
|----------|----------|
|----------|----------|


Helping all people  
live healthy lives

|                                         |                                                                                             |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                             |                                                                                             |
| Form:                                   | Liquid                                                                                      |
| Colour:                                 | Clear                                                                                       |
|                                         | Light red                                                                                   |
| Odour:                                  | Characteristic                                                                              |
| pH-value:                               | 7,5                                                                                         |
| Boiling point/Boiling range             | 82 °C                                                                                       |
| Flash point                             | 25 °C                                                                                       |
| Ignition temperature                    | 410,0 °C                                                                                    |
| Self-igniting                           | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 2,0 Vol %                                                                                   |
| Upper:                                  | 44,0 Vol %                                                                                  |
| Vapour pressure at 20 °C:               | 128,0 hPa                                                                                   |
| Solubility in / Miscibility with Water: | Soluble                                                                                     |

7. BD CytoRich Blue Preservative - plastic container, 3,6 L. High Density Polyethylene (HDPE) bottle and 38/400 pp White r/m f217 cap, proprietary formula for the preservative fluid CytoRich blue. It is an alcohol-based preservative fluid intended for preservation of cytology cells in suspension and their subsequent preparation for cytological evaluation.

#### Composition of BD CytoRich Blue Preservative

|            |         |       |
|------------|---------|-------|
| S: 64-17-5 | ethanol | 44,0% |
|------------|---------|-------|

Helping all people  
live healthy lives

|                                   |                                                                                                                                                  |      |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|------|
| EINECS: 200-578-6                 | F R10<br>Flam. Liq. 2, H225                                                                                                                      |      |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 5,0% |
| CAS: 25322-68-3                   | Polyethylene Glycol<br>Flam. Liq. 3, H226                                                                                                        | 1,0% |

### Physical and chemical properties of BD CytoRich Blue Preservative

| Property                                   | Comments                                                                                    |
|--------------------------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                                |                                                                                             |
| Form:                                      | Liquid                                                                                      |
| Colour:                                    | Clear<br>Blue                                                                               |
| Odour:                                     | Alcohol-like                                                                                |
| Boiling point/Boiling range:               | 82 °C                                                                                       |
| Flash point:                               | 15 °C                                                                                       |
| Ignition temperature:                      | 425,0 °C                                                                                    |
| Self-igniting:                             | Product is not selfigniting.                                                                |
| Danger of explosion                        | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                          |                                                                                             |
| Lower:                                     | 3,5 Vol %                                                                                   |
| Upper:                                     | 15,0 Vol %                                                                                  |
| Density:                                   | 0,93 g/cm <sup>3</sup>                                                                      |
| Solubility in / Miscibility with<br>Water: | Soluble                                                                                     |
| Solvent content:                           |                                                                                             |
| Organic solvents:                          | 49,0 %                                                                                      |
| Water:                                     | 49,0 %                                                                                      |

8. BD PrepStain Non-GYN Stain kit consists of Hematoxylin Stain 0,5 and EA/OG Combo Stain. Both stains are standard stains in cytology. High Density Polyethylene (HDPE) bottles and caps, proprietary formula for the Hematoxylin and EA/OG combo stain.

- Hematoxylin Stain 0,5 - 1 vial, 480 mL

### Physical and chemical properties of Hematoxylin Stain 0,5

| Property                                   | Comments                                      |
|--------------------------------------------|-----------------------------------------------|
| Appearance:                                |                                               |
| Form:                                      | Liquid                                        |
| Colour:                                    | Violet<br>Red                                 |
| Odour:                                     | Characteristic                                |
| pH-value:                                  | 2,6                                           |
| Flash point                                | >100 °C                                       |
| Ignition temperature:                      | 410,0 °C                                      |
| Self-igniting:                             | Product is not selfigniting.                  |
| Danger of explosion                        | Product does not present an explosion hazard. |
| Density:                                   | 1,02 g/cm <sup>3</sup>                        |
| Solubility in / Miscibility with<br>Water: | Soluble                                       |

- EA/OG Combo Stain - 1 vial, 480 mL

### Composition of EA/OG Combo Stain

|                                   |                                                                                                                                                  |       |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 64-17-5<br>EINECS: 200-578-6 | ethanol<br>F R10<br>Flam. Liq. 2, H225                                                                                                           | 69,7% |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 25,0% |
| CAS: 64-19-7<br>EINECS: 200-580-7 | acetic acid<br>C R35<br>R10<br>Flam. Liq. 3, H226; Skin Corr. 1A, H314                                                                           | 2,0%  |

### Physical and chemical properties of EA/OG Combo Stain

| Property                    | Comments                                                                                    |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                 |                                                                                             |
| Form:                       | Liquid                                                                                      |
| Colour:                     | Dark green                                                                                  |
| Odour:                      | Alcohol-like                                                                                |
| pH-value                    | 5,1 +/- 0,2                                                                                 |
| Boiling point/Boiling range | 73 °C                                                                                       |
| Flash point                 | 12 °C                                                                                       |
| Ignition temperature        | 425 °C                                                                                      |
| Self-igniting               | Product is not selfigniting.                                                                |
| Danger of explosion         | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:           |                                                                                             |
| Lower:                      | 3,5 Vol %                                                                                   |
| Upper:                      | 44,0 Vol %                                                                                  |
| Density:                    | 0,826 g/cm <sup>3</sup>                                                                     |

Helping all people  
live healthy lives

|                                         |                                  |
|-----------------------------------------|----------------------------------|
| Solubility in / Miscibility with Water: | Not miscible or difficult to mix |
| Solvent content:                        |                                  |
| Organic solvents:                       | 96,7 %                           |
| Water:                                  | 2,5 %                            |

9. BD SurePath Preservative Fluid in vials 10 mL (kit - 25 vials) is an alcohol-based, preservation solution that serves as a transport, preservative and antibacterial medium for gynecologic specimens.

#### Composition of BD SurePath Preservative Fluid

|                                   |                                                                                                                                                  |       |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 64-17-5<br>EINECS: 200-578-6 | ethanol<br>F R10<br>Flam. Liq. 2, H225                                                                                                           | 21,7% |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 1,2%  |
| CAS: 67-63-0<br>EINECS: 200-661-7 | isopropanol<br>Xi R36; F R11<br>R67<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                                   | 1,1%  |

#### Physical and chemical properties of BD SurePath Preservative Fluid


| Property | Comments |
|----------|----------|
|----------|----------|

Helping all people  
live healthy lives

|                                         |                                                                                             |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                             |                                                                                             |
| Form:                                   | Liquid                                                                                      |
| Colour:                                 | Clear                                                                                       |
|                                         | Colourless                                                                                  |
| Odour:                                  | Alcohol-like                                                                                |
| pH-value:                               | 7,2                                                                                         |
| Boiling point/Boiling range:            | 86 °C                                                                                       |
| Flash point                             | 35 °C                                                                                       |
| Ignition temperature:                   | 425,0 °C                                                                                    |
| Self-igniting                           | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 3,5 Vol %                                                                                   |
| Upper:                                  | 15,0 Vol %                                                                                  |
| Vapour pressure at 20 °C:               | 59,0 hPa                                                                                    |
| Solubility in / Miscibility with Water: | Soluble                                                                                     |
| Solvent content:                        |                                                                                             |
| Organic solvents:                       | 24,2 %                                                                                      |
| Water:                                  | 74,7 %                                                                                      |
| Solids content:                         | 0,9 %                                                                                       |

10. BD SurePath Preservative Fluid in vials 10 mL (kit - 500 vials) is an alcohol-based, preservation solution that serves as a transport, preservative and antibacterial medium for gynecologic specimens.

#### Composition of BD SurePath Preservative Fluid

|                   |                    |       |
|-------------------|--------------------|-------|
| CAS: 64-17-5      | ethanol            | 21,7% |
| EINECS: 200-578-6 | F R10              |       |
|                   | Flam. Liq. 2, H225 |       |

Helping all people  
live healthy lives

|                                   |                                                                                                                                                  |      |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|------|
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 1,2% |
| CAS: 67-63-0<br>EINECS: 200-661-7 | isopropanol<br>Xi R36; F R11<br>R67<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                                   | 1,1% |

#### Physical and chemical properties of BD SurePath Preservative Fluid

| Property                                | Comments                                                                                    |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Appearance:                             |                                                                                             |
| Form:                                   | Liquid                                                                                      |
| Colour:                                 | Clear<br>Colourless                                                                         |
| Odour:                                  | Alcohol-like                                                                                |
| pH-value:                               | 7,2                                                                                         |
| Boiling point/Boiling range:            | 86 °C                                                                                       |
| Flash point                             | 35 °C                                                                                       |
| Ignition temperature:                   | 425,0 °C                                                                                    |
| Self-igniting                           | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 3,5 Vol %                                                                                   |
| Upper:                                  | 15,0 Vol %                                                                                  |
| Vapour pressure at 20 °C:               | 59,0 hPa                                                                                    |
| Solubility in / Miscibility with Water: | Soluble                                                                                     |
| Solvent content:                        |                                                                                             |
| Organic solvents:                       | 24,2 %                                                                                      |
| Water:                                  | 74,7 %                                                                                      |
| Solids content:                         | 0,9 %                                                                                       |

11. BD SureDetect Counterstains consists of Hematoxylin Stain and Blue Stain. Plastic vials and caps, proprietary formula for the reagents. Hematoxylin Stain is Aqueous based Mayer's Hematoxylin. Blue Stain is Tris buffered saline, pH 7.4 with Tween 20 and 0.09% NaN<sub>3</sub>.

- Hematoxylin Stain in a vial 18 mL

Helping all people  
live healthy lives

### Composition of Hematoxylin Stain

|                   |                                                             |      |
|-------------------|-------------------------------------------------------------|------|
| CAS: 302-17-0     | chloral hydrate                                             | 5,0% |
| EINECS: 206-117-5 | T R25; Xi R36/38                                            |      |
|                   | Acute Tox. 3, H301; Skin Irrit. 2, H315; Eye Irrit. 2, H319 |      |

### Physical and chemical properties of Hematoxylin Stain

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Violet<br>Red                                 |
| Odour:                                  | Characteristic                                |
| pH-value                                | 2,4                                           |
| Self-igniting                           | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Vapour pressure:                        | 23,0 hPa                                      |
| Solubility in / Miscibility with Water: | Soluble                                       |

- Blue Stain in a vial 18 mL

### Physical and chemical properties of Blue Stain

| Property                                   | Comments                                      |
|--------------------------------------------|-----------------------------------------------|
| Appearance:                                |                                               |
| Form:                                      | Liquid                                        |
| Odour:                                     | Characteristic                                |
| Self-igniting                              | Product is not selfigniting.                  |
| Danger of explosion                        | Product does not present an explosion hazard. |
| Vapour pressure at 20 °C:                  | 23 hPa                                        |
| Density at 20 °C:                          | 1 g/cm <sup>3</sup>                           |
| Solubility in / Miscibility with<br>Water: | Soluble                                       |
| Solvent content:<br>Water:                 | 98,1 %                                        |

12. BD ProEx C Immunocytochemical Test consists of primary antibodies to protein supported MCM2 and to TOP2 $\alpha$ . Ab Cocktail is Ready-to-use monoclonal antibody cocktail supplied in TRIS buffered solution with Tween 20, pH 7.4. Contains stabilizing proteins and anti-microbial agent.

- vial 15 mL – primary antibodies to protein supported MCM2 and to TOP2 $\alpha$

Helping all people  
live healthy lives

**Physical and chemical properties of primary antibodies to protein supported  
MCM2 and to TOP2 $\alpha$**

| Property                                   | Comments                                      |
|--------------------------------------------|-----------------------------------------------|
| Appearance:                                |                                               |
| Form:                                      | Liquid                                        |
| Odour:                                     | Characteristic                                |
| pH-value:                                  | 7,4                                           |
| Self-igniting:                             | Product is not selfigniting.                  |
| Danger of explosion                        | Product does not present an explosion hazard. |
| Solubility in / Miscibility with<br>Water: | Soluble                                       |


13. BD SureDetect Detection Reagents are Peroxidase Block reagent, Protein Block reagent, Secondary Mouse Probe antibodies, Reagent polymer based on peroxidase, DAB Substrate Buffer (solvent for 3,3'-diaminobenzidine), DAB Chromagen (dye, 3,3'-diaminobenzidine). Plastic vials and caps, polyethylene pipette and graduated cylinder, proprietary formula for the reagents. Peroxidase Block reagent is buffered hydrogen peroxide plus stabilizer and proprietary components. Protein Block reagent is purified casein plus proprietary combination of proteins in modified PBS with preservative and surfactant. Secondary Mouse Probe antibodies are binds to mouse monoclonal antibodies. Reagent polymer based on peroxidase is Polymer conjugated with horseradish peroxidase that binds to Mouse Probe Reagent. DAB Substrate Buffer (solvent for 3,3'-diaminobenzidine) is substrate buffer used in the preparation of the DAB chromogen. DAB Chromagen is 3,3'-diaminobenzidine chromogen solution.

- vial 15 mL - Peroxidase Block reagent


Helping all people  
live healthy lives

#### Physical and chemical properties of Peroxidase Block reagent

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Clear                                         |
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7,0                                           |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |


- vial 15 mL – Protein Block reagent

### Physical and chemical properties of Protein Block reagent

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Blue                                          |
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7,0                                           |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

- vial 15 mL – Secondary Mouse Probe antibodies

### Physical and chemical properties of Secondary Mouse Probe antibodies

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Green                                         |
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7,0                                           |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

- vial 15 mL – Reagent polymer based on peroxidase


Helping all people  
live healthy lives

### Physical and chemical properties of Reagent polymer based on peroxidase

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Light brown<br>Green                          |
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7.0                                           |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

- vial 18 mL – DAB Substrate Buffer (solvent for 3,3'- diaminobenzidine)

### Composition of DAB Substrate Buffer (solvent for 3,3'- diaminobenzidine)

Helping all people  
live healthy lives

|                                   |                                                                                                     |      |
|-----------------------------------|-----------------------------------------------------------------------------------------------------|------|
| CAS: 91-95-2<br>EINECS: 202-110-6 | 3,3 diaminobenzidine<br>T Carc. Cat. 2 R45; Xn R68<br>Muta. Cat. 3<br>Muta. 2, H341; Carc. 1B, H350 | 1,0% |
|-----------------------------------|-----------------------------------------------------------------------------------------------------|------|

### Physical and chemical properties of DAB Substrate Buffer (solvent for 3,3'-diaminobenzidine)

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Colour:                                 | Clear                                         |
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7,0                                           |
| Self-igniting:                          | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

- vial 1.0 mL - DAB Chromagen (dye, 3,3'- diaminobenzidine)

### Composition of DAB Chromagen (dye, 3,3'- diaminobenzidine)

|                                   |                                                                                                     |      |
|-----------------------------------|-----------------------------------------------------------------------------------------------------|------|
| CAS: 91-95-2<br>EINECS: 202-110-6 | 3,3 diaminobenzidine<br>T Carc. Cat. 2 R45; Xn R68<br>Muta. Cat. 3<br>Muta. 2, H341; Carc. 1B, H350 | 1,0% |
|-----------------------------------|-----------------------------------------------------------------------------------------------------|------|

### Physical and chemical properties of DAB Chromagen (dye, 3,3'-diaminobenzidine)

| Property    | Comments |
|-------------|----------|
| Appearance: |          |
| Form:       | Liquid   |
| Colour:     | Brown    |

Helping all people  
live healthy lives

|                                         |                                               |
|-----------------------------------------|-----------------------------------------------|
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7,0                                           |
| Self-igniting:                          | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

14. BD SureDetect SiHa Cell Control in a vial 8 mL for preparation of samples. Low Density Polyethylene (LDPE) bottle and tip, proprietary formula for the reagents.

#### Composition of BD SureDetect SiHa Cell Control

|                                   |                                                                                                                                                  |       |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| CAS: 64-17-5<br>EINECS: 200-578-6 | ethanol<br>F R10<br>Flam. Liq. 2, H225                                                                                                           | 21,7% |
| CAS: 67-56-1<br>EINECS: 200-659-6 | methanol<br>T R23/24/25-39/23/24/25; F R11<br>Flam. Liq. 2, H225; Acute Tox. 3, H301; Acute Tox. 3, H311;<br>Acute Tox. 3, H331; STOT SE 1, H370 | 1,2%  |
| CAS: 67-63-0<br>EINECS: 200-661-7 | isopropanol<br>Xi R36; F R11<br>R67<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                                   | 1,1%  |

#### Physical and chemical properties of BD SureDetect SiHa Cell Control


| Property    | Comments     |
|-------------|--------------|
| Appearance: |              |
| Form:       | Liquid       |
| Colour:     | Clear        |
| Odour:      | Alcohol-like |
| pH-value:   | 7,2          |

Helping all people  
live healthy lives

|                                         |                                                                                             |
|-----------------------------------------|---------------------------------------------------------------------------------------------|
| Boiling point/Boiling range:            | 86 °C                                                                                       |
| Flash point:                            | 35 °C                                                                                       |
| Ignition temperature:                   | 425,0 °C                                                                                    |
| Self-igniting:                          | Product is not selfigniting.                                                                |
| Danger of explosion                     | Product is not explosive. However, formation of explosive air/vapour mixtures are possible. |
| Explosion limits:                       |                                                                                             |
| Lower:                                  | 3,5 Vol %                                                                                   |
| Upper:                                  | 15,0 Vol %                                                                                  |
| Vapour pressure:                        | 59,0 hPa                                                                                    |
| Solubility in / Miscibility with Water: | Soluble                                                                                     |
| Solvent content:                        |                                                                                             |
| Organic solvents:                       | 24,2 %                                                                                      |

15. BD SureDetect Slide Preparation Buffer (10x). High Density Polyethylene (HDPE) bottle and cap, proprietary formula for the reagents. BD SureDetect Slide Preparation Buffer is to be used on cytology and histology specimens, which have been mounted on glass for target retrieval prior to immunocytochemistry (ICC) or immunohistochemistry (IHC) procedures.

#### Physical and chemical properties of BD SureDetect Slide Preparation Buffer

| Property       | Comments                     |
|----------------|------------------------------|
| Appearance:    |                              |
| Form:          | Liquid                       |
| Colour:        | Clear                        |
| Odour:         | Characteristic               |
| pH-value:      | 6,5                          |
| Self-igniting: | Product is not selfigniting. |


Helping all people  
live healthy lives

|                                         |                                               |
|-----------------------------------------|-----------------------------------------------|
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

16. BD ProEx C Immunohistochemical Test in a vial 7 mL contains an antibody reagent for use in a three-step manual or automated immunohistochemical staining procedure in formalin-fixed paraffin-embedded tissue biopsy specimens.

#### Physical and chemical properties of BD ProEx Immunohistochemical Test

| Property                                | Comments                                      |
|-----------------------------------------|-----------------------------------------------|
| Appearance:                             |                                               |
| Form:                                   | Liquid                                        |
| Odour:                                  | Characteristic                                |
| pH-value:                               | 7,4                                           |
| Self-igniting:                          | Product is not selfigniting.                  |
| Danger of explosion                     | Product does not present an explosion hazard. |
| Solubility in / Miscibility with Water: | Soluble                                       |

#### 8. Safety precautions

- For in vitro Diagnostic Use.
- Cytologic specimens may contain infectious agents. Follow appropriate biohazard precautions when handling samples.
- For professional use only.
- Good laboratory practices should be followed.
- Avoid splashing or generating aerosols. Operators should use appropriate hand, eye/face and clothing protection.
- Do not pipette by mouth.
- Do not ingest (contains denatured alcohol).

- 3,3'-Diaminobenzidine (DAB) is classified as a suspected carcinogen with the potential risk of producing irreversible effects. Avoid physical contact and prolonged or repeated exposure. A chemical fume hood is highly recommended. Do not get in eyes, on skin or on clothing. Wash hands thoroughly after handling.
- Some reagents contain sodium azide ( $\text{NaN}_3$ ), a highly toxic chemical in pure form. Sodium azide may be fatal if swallowed or absorbed through skin. These reagents are harmful if inhaled. These reagents are hazardous solids. Causes irritation to skin, eyes and respiratory tract. Affects central nervous system, kidneys and cardiovascular system. Although not classified as hazardous at product concentrations, build-ups of sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing.
- Hematoxylin and Blue reagent: These reagents are harmful if swallowed. Hematoxylin and Blue reagents are eye, skin and respiratory irritants.
- Minimize microbial contamination of reagents to avoid nonspecific staining.
- Incubation times, temperatures or methods other than those specified may give erroneous results.
- Reagents have been diluted for optimal performance. Further dilution may result in loss of antigen staining.

## 9. Packaging

Bottles of reagents with instructions for use are placed in a box. Material of the box is cardboard.

## 10. Labeling

| Symbol                                                                              | Description                        |
|-------------------------------------------------------------------------------------|------------------------------------|
|  | Manufacturer                       |
|  | European representative            |
|  | Code number                        |
|  | Serial number                      |
|  | CE marking                         |
|  | Read the instructions before use   |
| IVD                                                                                 | In Vitro Diagnostic Medical Device |

Helping all people  
live healthy lives

|                                                                                     |                                                  |
|-------------------------------------------------------------------------------------|--------------------------------------------------|
|    | Temperature limitation                           |
|    | Use by                                           |
|    | Warning                                          |
|    | Danger. May cause damage to organs.              |
|   | Danger. Highly Flammable liquid and vapour.      |
|  | Danger. Causes severe skin burns and eye damage. |
|  | Contains sufficient for <n> tests                |
|  | Single use                                       |
|  | Fragile, handle with care                        |

### 11. Transportation, storage and shelf life

The transportation is carried out by all types of covered transport at temperatures from 15° C to 30° C and relative humidity of no more than 90%, non-condensing, in accordance with the shipping rules applicable to the kind of transport with the exception of BD SureDetect Counterstains, BD ProEx C Immunohistochemical Test, BD SureDetect Detection Reagents and BD SureDetect SiHa Cell Control for which transportation is carried out at temperatures from 2° C to 8° C.

Store in cool, dry conditions in well sealed receptacles; protect from direct sunlight. Do not store together with oxidizing and acidic materials. Do not store together with

alkalis (caustic solutions).

Shelf life is only applicable for:

Alcohol Blend Rinse – 36 months

Hematoxylin Stain 0.5 – 12 months

EA/OG Combo Stain - 12 months

BD SurePath Preservative Fluid – 36 months

Hematoxylin Stain 0.75 – 12 months

Primary antibodies to protein supported MCM2 and to TOP2 $\alpha$  – 12 months

SiHa Cell Control – 9 months

Slide Preparation Buffer – 9 months

## 12. Sterilization

Reagents for preparation and coloring samples for cytological tests are not sterile.

## 13. Technical service maintenance and repair

Technical service maintenance and repair is not applicable to reagents for preparation and coloring samples for cytological tests.

## 14. Requirements for environmental protection during the use of the devices

Reagents for preparation and coloring samples for cytological tests in the process of the life cycle have no negative impact on the environment.

## 15. Disposal rules

Disposal of reagents for preparation and coloring samples for cytological tests must be in accordance with the laws and local regulations, as well as with the protocols of the medical institution.

It is recommended to recycle the reagents in compliance with Public health regulations (SanPiN) 2.1.7.2790-10 "Sanitary and epidemiological requirements to medical wastes treatment".

## 16. Manufacturer's warranty

Proper maintenance is necessary for the reagents for cytological tests. TriPath Imaging is not responsible, under warranty or otherwise, for damage due to abuse or neglect.

Helping all people  
live healthy lives

### 17. Denomination and legal address of the authorized representative in Russia

"Becton Dickinson Vostok" Limited Liability Company.

("Becton Dickinson Vostok" LLC)

Legal address: Block C, bldg. 1, 12, ul. Dvintsev, Moscow, Russian Federation, 127018

Actual address: Block C, bldg. 1, 12, ul. Dvintsev, Moscow, Russian Federation, 127018

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М.П. / Stamp

|                                          |                               |
|------------------------------------------|-------------------------------|
| State of                                 | <u>NORTH CAROLINA</u>         |
| County of                                | <u>DURHAM</u>                 |
| Sworn to and subscribed before me on the | <u>23</u> day of <u>MARCH</u> |
| Notary Public's Signature                | <u>[Signature]</u>            |
| My Commission Expires                    | <u>12/28/2018</u>             |

Helping all people  
live healthy lives

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