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Stamp

Besmed 

Instructions for use of the medical device:
Breathing anaesthetic circuit in design variants

1. Medical device name

Breathing anaesthetic circuit

Design variants:

- I. Breathing anaesthetic circuit ribbed for adults CT-16350, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Y-shaped connector – 1 piece.
 4. Tube outlet and inlet valve with or without Luer-lock port – 1 piece.
 5. Instructions for use – 1 piece.
- II. Breathing anaesthetic circuit ribbed for adults CT-16352, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Y-shaped connector – 1 piece.
 4. Tube outlet and inlet valve with or without Luer-lock port – 1 piece.
 5. Breathing latex bag – 1 piece.
 6. Instructions for use – 1 piece.
- III. Breathing anaesthetic circuit ribbed for adults CT-16370, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Y-shaped connector – 1 piece.
 4. Tube outlet and inlet valve with or without Luer-lock port – 1 piece.
 5. Instructions for use – 1 piece.
- IV. Breathing anaesthetic circuit ribbed for adults CT-16306, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Y-shaped connector – 1 piece.
 4. Tube outlet and inlet valve with or without Luer-lock port – 1 piece.
 5. Breathing bag – 1 piece.
 6. Instructions for use – 1 piece.
- V. Breathing anaesthetic circuit ribbed with a water collector for adults CT-16130, consisting of:
 1. Tube- 3 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Water collector – 1 piece.
 5. Connector 22M/22F – 4 pcs.
 6. Straight connector 22M/15F – 1 piece.
 7. Hook – 1 piece.
 8. Instructions for use – 1 piece.
- VI. Breathing anaesthetic circuit ribbed with a water collector for adults CT-16390, consisting of:
 1. Tube- 5 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 4. Water collector – 2 pcs.
 5. Connector 22M/22F – 4 pcs.
 6. Straight connector 22M/15F – 1 piece.
 7. Hook – 1 piece.
 8. Instructions for use – 1 piece.
- VII. Breathing anaesthetic circuit ribbed for children CT-18601, consisting of:
 1. Tube- 1 piece
 2. Tube outlet and inlet valve with Luer-lock port – 1 piece.

3. Instructions for use - 1 piece.
- VIII. Breathing anaesthetic circuit ribbed for children CT-18604, consisting of:
 1. Tube- 2 pcs
 2. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 3. Breathing bag – 1 piece.
 4. Straight connector 15F/22M – 1 piece.
 5. Instructions for use - 1 piece.
- IX. Breathing anaesthetic circuit riffled for adults CT-11335, consisting of:
 1. Tube- 2 pcs.
 2. Connector 22M/22F – 1 piece.
 3. Tube clamp – 2 pcs.
 4. Nebulizer – 1 piece.
 5. Heated wire – 1 piece.
 6. Instructions for use - 1 piece.
- X. Breathing anaesthetic circuit stretchable for adults CT-16012, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 4. Instructions for use - 1 piece.
- XI. Breathing anaesthetic circuit stretchable for adults CT-16014, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 4. Tube 1.5m – 1 piece.
 5. Breathing latex bag – 1 piece.
 6. Straight connector 15F/22M – 1 piece.
 7. Instructions for use - 1 piece.
- XII. Breathing anaesthetic circuit smoothbore for adults CT-10877, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Instructions for use - 1 piece.
- XIII. Breathing anaesthetic circuit smoothbore for adults CT-10867, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Breathing latex bag – 1 piece.
 5. Tube 1 m – 1 piece.
 6. Instructions for use - 1 piece.
- XIV. Breathing anaesthetic circuit smoothbore for adults CT-10879, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Instructions for use - 1 piece.
- XV. Breathing anaesthetic circuit smoothbore for adults CT-10869, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Breathing latex bag – 1 piece.
 5. Tube 1 m – 1 piece.
 6. Instructions for use - 1 piece.
- XVI. Breathing anaesthetic circuit smoothbore with water collector for adults CT-10889, consisting of:
 1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.

3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Water collector – 1 piece.
 5. Instructions for use - 1 piece.
- XVII. Breathing anaesthetic circuit smoothbore with water collector for adults CT-10899, consisting of:**
1. Tube- 2 pcs.
 2. Parallel connector – 1 piece.
 3. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 4. Water collector – 2 piece.
 5. Instructions for use - 1 piece.
- XVIII. Breathing anaesthetic circuit coaxial for adults CT-11110, consisting of:**
1. Tube- 1 piece.
 2. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 3. Instructions for use - 1 piece.
- XIX. Breathing anaesthetic circuit coaxial for adults CT-11170, consisting of:**
1. Tube- 1 piece.
 2. Tube outlet and inlet valve without Luer-lock port – 1 piece.
 3. CO2 tube – 1 piece.
 4. Instructions for use - 1 piece.
- XX. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19021, consisting of:**
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Instructions for use - 1 piece.
- XXI. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19022, consisting of:**
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Water collector – 1 piece.
 7. Instructions for use - 1 piece.
- XXII. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19023, consisting of:**
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Automatic feed chamber – 1 piece.
 7. Instructions for use - 1 piece.
- XXIII. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19024, consisting of:**
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Automatic feed chamber – 1 piece.
 7. Water collector – 1 piece.
 8. Instructions for use - 1 piece.
- XXIV. Breathing anaesthetic circuit with heated breathing mixture for children CT-19321, consisting of:**
1. Tube- 1 piece.

2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Water collector – 1 piece.
 7. Instructions for use - 1 piece.
- XXV. Breathing anaesthetic circuit with heated breathing mixture for children CT-19322, consisting of:
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Instructions for use - 1 piece.
- XXVI. Breathing anaesthetic circuit with heated breathing mixture for children CT-19323, consisting of:
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Automatic feed chamber – 1 piece.
 7. Water collector – 1 piece.
 8. Instructions for use - 1 piece.
- XXVII. Breathing anaesthetic circuit with heated breathing mixture for children CT-19324, consisting of:
1. Tube- 1 piece.
 2. Heated wire – 1 piece.
 3. Connector 15F – 1 piece.
 4. Heated inhalation tube – 1 piece.
 5. Heated exhalation tube – 1 piece.
 6. Automatic feed chamber – 1 piece.
 7. Instructions for use - 1 piece.
- XXVIII. Breathing anaesthetic circuit with heated breathing mixture for newborns CT-19240, consisting of:
1. Tube- 1 piece.
 2. Connector 10M – 1 piece.
 3. Heated inhalation tube – 1 piece.
 4. Heated exhalation tube – 1 piece.
 5. Instructions for use - 1 piece.
- XXIX. Breathing anaesthetic circuit with heated breathing mixture for newborns CT-19340, consisting of:
1. Tube- 1 piece.
 2. Connector 15F – 1 piece.
 3. Heated inhalation tube – 1 piece.
 4. Heated exhalation tube – 1 piece.
 5. Water collector – 1 piece.
 6. Instructions for use - 1 piece.
- XXX. Breathing anaesthetic circuit for adults CT-11730, consisting of:
1. Tube- 1 piece.
 2. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 3. Instructions for use - 1 piece.
- XXXI. Breathing anaesthetic circuit for adults CT-11731, consisting of:
1. Tube- 1 piece.
 2. Breathing bag – 1 piece.
 3. Straight connector – 1 piece.
 4. Tube outlet and inlet valve with Luer-lock port – 1 piece.
 5. Tube 90 cm – 1 piece.
 6. Instructions for use - 1 piece.

2. Purpose of the medical device

Breathing circuits for anesthesia (hereinafter referred to as circuits or breathing circuits) are intended for inhalation anesthesia. Inhalation anesthesia is a method of general anesthesia based on the use of gaseous or volatile general anesthetics that enter the patient's body through the respiratory tract. Breathing circuit – BDC is a set of devices that ensure the delivery of anaesthesia-respiratory mixture to the patient.

3. Scope of the medical device

in medical and health care facilities.

4. Medical device usage conditions

Temperature range: 10 ~ 40°C

Relative humidity range: 20% to 95%, non-condensing

Ambient pressure range: from 500 hPa to 1200 hPa

Background sound pressure level: <70 dBa

5. Indications for use

- for large surgical operations.
- long-term procedures and operations, the need preservation of anaesthetic, the ability to monitor ECG and saturation.

6. Contraindications/possible side effects

Absolute contraindications include:

- pathology of the cardiovascular system in the decompensation stage, malignant hypertension, severe hypotension caused by shock, collapse, intoxication, anaemia;
- pathology of the respiratory system with pronounced respiratory failure – acute pneumonia, bronchiectatic diseases, severe forms of tuberculosis;
- pathology of the hepatobiliary system with acute liver dysfunction, development of liver failure – cirrhosis, hepatitis;
- kidney pathology with severe impairment of their functions, development of renal failure – severe glomerulonephritis;
- severe endocrine pathologies – diabetes mellitus in the stage of decompensation of metabolic processes;
- severe anaemias;
- severe form of cachexia;
- increased intracranial pressure in brain tumors, hematomas.

Relative contraindications are various pathologies of organs and systems in which violations of their functions are expressed slightly and are subject to medication correction.

7. Warnings

- Do not use for transfusion of blood or blood components.
- Do not use if the package is opened or damaged.
- In the absence of an alternative method of respiratory support for the patient and the necessary spare equipment.
- In the absence of an adequate source of compressed medical gases.
- When transporting patients outside the clinic.
- due to the presence of valves in the circuit, patients with low body weight and/or reduced strength of respiratory movements may have difficulty overcoming the resistance of the circuit, which leads to hypoventilation.

8. Safety instructions


- The procedure should be performed according to the attending physician's prescription.
 - Intended for single use only. Re-use is prohibited.
 - When preparing for work, check the breathing circuit for the presence of foreign bodies.
 - The procedure should be performed by qualified medical personnel only.
 - For use only for one patient
- Inadequate anesthesia may lead to coughing, breath retention, or laryngeal spasm.

9. Medical device manufacturer

Besmed Health Business Corp
(No. 5, Lane 116, Wu-Kong 2nd Rd,
Wu-Ku District, New Taipei City 248, Taiwan)
Tel.: 886-2-2290-3959
Fax: 886-2-2299-9076
E-mail address: info@besmed.com
Web-site: www.besmed.com

10. Basic composition (completeness) of the medical device, the list of design variants

Breathing systems, ribbed for adults

SPEC.


Patient	15F/22M W/ Luer
Elbow	22M

Type	Length	1	2	3	4	CODE
Parallel 	1.5m		•			CT-16350
						CT-16352
				•		CT-16354
		•	•		•	CT-16360
		•		•	•	CT-16362
	1.8m		•			CT-16370
						CT-16372
				•		CT-16306
		•	•		•	CT-16378
		•		•	•	CT-16380


Model	Description
Breathing anaesthetic circuit ribbed for adults CT-16350	Adult parallel breathing circuit, with tube length 1.5 m, diameter 22 mm (2 pcs), outlet and inlet valve 15F/22M, parallel connector 22M, Y-shaped connector 22M
Breathing anaesthetic circuit ribbed for adults CT-16352	Adult parallel breathing circuit, with tube length 1.5 m, diameter 22 mm (2 pcs), outlet and inlet valve 15F/22M, parallel connector 22M, Y-shaped connector 22M, breathing latex bag 2 L
Breathing anaesthetic circuit ribbed for adults CT-16370	Adult parallel breathing circuit, with tube length 1.8 m, diameter 22 mm (2 pcs), outlet and inlet valve 15F/22M, parallel connector 22M, Y-shaped connector 22M
Breathing anaesthetic circuit ribbed for adults CT-16306	Adult parallel breathing circuit with tube diameter 22 mm, length 1.8 m (2 pcs), breathing bag 2 L, tube 1 m, parallel connector 22M, Y-shaped connector 22M, outlet and inlet valve 15F/22M
Breathing anaesthetic circuit ribbed with a water collector for adults CT-16130	tube 75 cm-2 pcs, tube 1.5 m, diameter 22 mm, straight connector 22M/15F, hook connector 22M/22F-4 pcs, outlet and inlet valve without Luer-lock port 15F/22M, water collector, parallel connector 22M
Breathing anaesthetic circuit ribbed with a water collector for adults CT-16390	parallel connector 22M, tube 75 cm-4 pcs, water collector-2 pcs, tube 45 cm, straight connector 22M/15F, hook, connectors 22M/22F-4 pcs, outlet and inlet valve with Luer-lock port 15F/22M

Ribbed breathing circuits for children

Breathing anaesthetic circuit ribbed for children CT-18601	Breathing circuit tube diameter 10 mm, length 1.8 m, tube inlet and outlet valve with Luer-lock port 15F/22M
--	--

Breathing anaesthetic circuit ribbed for children CT-18604	Breathing circuit tube diameter 10 mm, tube length 1.8 m, tube inlet and outlet valve with Luer-lock port 15F/22M, breathing bag 0.5 L, tube 1 m, straight connector 15F/22M
Riffled breathing circuit	
	
Breathing anaesthetic circuit riffled for adults CT-11335	Breathing circuit with tubes - 15 cm and 100 cm, nebulizer 20 cc, connector 22M/22F, heated wire 110 cm, tube clip - 2 pcs.
Stretching breathing circuit	
Breathing anaesthetic circuit stretchable for adults CT-16012	Adult breathing circuit, parallel, tube length 2.0 m, tube diameter 22 mm (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve with Luer-lock port 15F/22M
Breathing anaesthetic circuit stretchable for adults CT-16014	Adult breathing circuit, parallel tube diameter 22 mm, tube length 2.0 m (2 pcs), tube 1.5 m, parallel connector 22M/22M, tube inlet and outlet valve with Luer-lock port 15F/22M, straight connector 15F/22M, breathing latex bag 2 L
Smooth-bore breathing circuit	
Breathing anaesthetic circuit smooth-bore for adults CT-10877	Adult parallel breathing circuit, tube diameter 22 mm, tube length 1.5 m (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve without Luer-lock port 15F/22M
Breathing anaesthetic circuit smooth-bore for adults CT-10867	Adult parallel breathing circuit tube diameter 22 mm, tube length 1.5 m (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve without Luer-lock port 15F/22M, latex bag 2 L, tube 1 m
Breathing anaesthetic circuit smooth-bore for adults CT-10879	Adult parallel breathing circuit, tube diameter 22 mm, tube length 1.8 m (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve without Luer-lock port 15F/22M
Breathing anaesthetic circuit smooth-bore for adults CT-10869	Adult parallel breathing circuit tube diameter 22 mm, tube length 1.8 m (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve without Luer-lock port 15F/22M, latex bag 2 L, tube 1 m
Breathing anaesthetic circuit smooth-bore with moisture collector for adults CT-10889	Adult parallel breathing circuit, tube diameter 22 mm, tube length 1.5 m (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve without Luer-lock port 15F/22M, water collector
Breathing anaesthetic circuit smooth-bore with moisture collector for adults CT-10899	Adult parallel breathing circuit tube diameter 22 mm, tube length 1.5 m (2 pcs), parallel connector 22M/22M, tube inlet and outlet valve without Luer-lock port 15F/22M, water collector - 2 pcs

Coaxial breathing circuit


Breathing anaesthetic circuit coaxial for adults CT-11110

Coaxial circuit for anesthesia tube len 1.5 m, tube inlet and outlet valve with lock port 15F/22M

Breathing anaesthetic circuit coaxial for adults CT-11170

Coaxial circuit for anesthesia tube len 1.5 m, tube inlet and outlet valve with Luer-lock port 15F/22M, CO2 tube

Breathing circuit with heated breathing mixture for adults

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
Breathing anaesthetic circuit with heated breathing mixture for adults CT-19021

Breathing circuit for 2 connections tu diameter 22 mm, length 60 cm 22F/2 heated wire – 2 m, connector with Lu lock port 15F, heated inhalation tube

	22F/22F, heated exhalation tube 1.5 m + 30 cm 22F/22F
Breathing anaesthetic circuit with heated breathing mixture for adults CT-19022	Breathing circuit for 1 connection tube diameter 22 mm, length 60 cm 22F/22F, heated wire – 2m, connector with Luer-lock port 15F, heated inhalation tube 1.5 m 22F/22F, heated exhalation tube 30+120+175 cm 22F/22F, water collector
Breathing anaesthetic circuit with heated breathing mixture for adults CT-19023	Breathing circuit for 2 connections tube diameter 22 mm, length 60 cm 22F/22F, heated wire – 2 m, connector with Luer-lock port 15F, heated inhalation tube 1.5 m 22F/22F, heated exhalation tube 1.5 m + 30 cm 22F/22F, automatic feed chamber
Breathing anaesthetic circuit with heated breathing mixture for adults CT-19024	Breathing circuit for 1 connection tube diameter 22 mm, length 60 cm 22F/22F, heated wire – 2m, connector with Luer-lock port 15F, heated inhalation tube 1.5 m 22F/22F, heated exhalation tube 30+120+175 cm 22F/22F, water collector, automatic feed chamber
Breathing circuit with heated breathing mixture for children	
Breathing anaesthetic circuit with heated breathing mixture for children CT-19321	Breathing circuit for 1 connection tube diameter 10 mm, length 60 cm 22F/15M, heated wire – 1.8 m, connector 15F, inhalation tube 110 cm+ 35 cm 22F, heated exhalation tube 90 cm+60 cm 15M, water collector
Breathing anaesthetic circuit with heated breathing mixture for children CT-19322	Breathing circuit for 2 connections tube diameter 10 mm, length 60 cm 22F/15M, heated wire – 1.8 m, connector 15F, heated inhalation tube 110 cm + 35 cm 22F, heated exhalation tube 150 cm 15M
Breathing anaesthetic circuit with heated breathing mixture for children CT-19323	Breathing circuit for 1 connection tube diameter 10 mm, length 60 cm 22F/15M, heated wire – 1.8 m, connector 15F, heated inhalation tube 110 cm+ 35 cm 22F, heated exhalation tube 90 cm+60 cm 15M, water collector, automatic feed chamber
Breathing anaesthetic circuit with heated breathing mixture for children CT-19324	Breathing circuit for 2 connections tube diameter 10 mm, length 60 cm 22F/15M, heated wire – 1.8 m, connector 15F, heated inhalation tube 110 cm+ 35 cm 22F, heated exhalation tube 150 cm 15M, automatic feed chamber
Breathing circuit with heated breathing mixture for newborns	

	CODE CT-19240
	CODE CT-19340
Breathing anaesthetic circuit with heated breathing mixture for newborns CT-19240	Breathing circuit tube diameter 10 mm heated inhalation tube 110 cm 22F, he wire – 1.8 m, connector 10M/10M, he exhalation tube 150 cm 10F/15M
Breathing anaesthetic circuit with heated breathing mixture for newborns CT-19340	Breathing circuit tube diameter 10 mm 60 cm 22F/10F, heated inhalation tube 110+35 cm 10F/22F, connector 15F, h exhalation tube 90+60 cm 10F/10F, w collector
Breathing anaesthetic circuit for adults	
Breathing anaesthetic circuit for adults CT-11730	Breathing circuit tube diameter 22 mm length 2.0 m, tube inlet and outlet valve with Luer-lock port 15F/22M
Breathing anaesthetic circuit for adults CT-11731	Breathing circuit tube diameter 22 mm length 2.0 m, tube inlet and outlet valve with Luer-lock port 15F/22M, breathing 2 L, straight connector 22M/15F, tube cm

MD package contents:

1. Breathing anaesthetic circuit ribbed for adults CT-16350:
- tube- 2 pcs.

- parallel connector — 1 piece
- Y-shaped connector — 1 piece
- tube outlet and inlet valve with or without Luer-lock port — 1 piece

2. Breathing anaesthetic circuit ribbed for adults CT-16352

- tube- 2 pcs.
- parallel connector — 1 piece
- Y-shaped connector — 1 piece
- tube outlet and inlet valve with or without Luer-lock port — 1 piece
- breathing latex bag — 1 piece

3. Breathing anaesthetic circuit ribbed for adults CT-16370

- tube- 2 pcs.
- parallel connector — 1 piece
- Y-shaped connector — 1 piece
- tube outlet and inlet valve with or without Luer-lock port — 1 piece

4. Breathing anaesthetic circuit ribbed for adults CT-16306

- tube- 2 pcs.
- parallel connector — 1 piece
- Y-shaped connector — 1 piece
- tube outlet and inlet valve with or without Luer-lock port — 1 piece
- breathing bag — 1 piece

10 mm,)
2F, heat
OM, heat
M
10 mm, t
on tube
15F, heat
10F, water

5. Breathing anaesthetic circuit ribbed with a water collector for adults CT-16130

- tube- 3 pcs.
- parallel connector — 1 piece
- tube outlet and inlet valve without Luer-lock port — 1 piece
- water collector — 1 piece
- connector 22M/22F — 4 pcs.
- straight connector 22M/15F — 1 pcs.
- hook — 1 piece

6. Breathing anaesthetic circuit ribbed with a water collector for adults CT-16390

- tube- 5 pcs.
- parallel connector — 1 piece
- tube outlet and inlet valve with Luer-lock port — 1 piece
- water collector — 2 piece
- connector 22M/22F — 4 pcs.
- straight connector 22M/15F — 1 pcs.
- hook — 1 piece

7. Breathing anaesthetic circuit ribbed for children CT-18601

- tube- 1 piece
- tube outlet and inlet valve with Luer-lock port — 1 piece

22 mm, t
et valve

8. Breathing anaesthetic circuit ribbed for children CT-18604

- tube- 2 pcs.
- tube outlet and inlet valve with Luer-lock port — 1 piece
- breathing bag — 1 piece
- straight connector 15F/22M — 1 piece

22 mm, t
et valve
breathing
, tube 90

9. Breathing anaesthetic circuit riffled for adults CT-11335

- tube- 2 pcs.
- connector 22M/22F — 1 pcs.
- tube clamp — 2 pcs.

- nebulizer — 1 piece
 - heated wire — 1 piece
10. Breathing anaesthetic circuit stretchable for adults CT-16012
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve with Luer-lock port — 1 piece
11. Breathing anaesthetic circuit stretchable for adults CT-16014
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve with Luer-lock port — 1 piece
 - tube- 1.5 m.
 - breathing latex bag — 1 piece
 - straight connector 15F/22M — 1 piece
12. Breathing anaesthetic circuit smooth-bore for adults CT-10877
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve without Luer-lock port — 1 piece
13. Breathing anaesthetic circuit smooth-bore for adults CT-10867
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve without Luer-lock port — 1 piece
 - breathing latex bag — 1 piece
 - tube 1 m — 1 piece
14. Breathing anaesthetic circuit smooth-bore for adults CT-10879
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve without Luer-lock port — 1 piece
15. Breathing anaesthetic circuit smooth-bore for adults CT-10869
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve without Luer-lock port — 1 piece
 - breathing latex bag — 1 piece
 - tube 1 m — 1 piece
16. Breathing anaesthetic circuit smooth-bore with moisture collector for adults CT-10889
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve without Luer-lock port — 1 piece
 - water collector — 1 piece
17. Breathing anaesthetic circuit smooth-bore with moisture collector for adults CT-10899
- tube- 2 pcs.
 - parallel connector — 1 piece
 - tube outlet and inlet valve without Luer-lock port — 1 piece
 - water collector — 2 piece
18. Breathing anaesthetic circuit coaxial for adults CT-11110
- tube- 1 piece
 - tube outlet and inlet valve with Luer-lock port — 1 piece
19. Breathing anaesthetic circuit coaxial for adults CT-11170

- tube- 1 piece
- tube outlet and inlet valve without Luer-lock port — 1 piece
- CO2 tube — 1 piece

20. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19021

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece

21. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19022

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- water collector — 1 piece

22. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19023

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- automatic feed chamber — 1 piece

23. Breathing anaesthetic circuit with heated breathing mixture for adults CT-19024

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- automatic feed chamber — 1 piece
- water collector — 1 piece

24. Breathing anaesthetic circuit with heated breathing mixture for children CT-19321

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- water collector — 1 piece

25. Breathing anaesthetic circuit with heated breathing mixture for children CT-19322

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece

26. Breathing anaesthetic circuit with heated breathing mixture for children CT-19323

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece

- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- automatic feed chamber — 1 piece
- water collector — 1 piece

27. Breathing anaesthetic circuit with heated breathing mixture for children CT-19324

- tube- 1 piece
- heated wire — 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- automatic feed chamber — 1 piece

28. Breathing anaesthetic circuit with heated breathing mixture for newborns CT-19240

- tube- 1 piece
- connector 10M — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece

29. Breathing anaesthetic circuit with heated breathing mixture for newborns CT-19340

- tube- 1 piece
- connector 15F — 1 piece
- heated inhalation tube — 1 piece
- heated exhalation tube — 1 piece
- water collector — 1 piece

30. Breathing anaesthetic circuit for adults CT-11730

- tube- 1 piece
- tube outlet and inlet valve with Luer-lock port — 1 piece

31. Breathing anaesthetic circuit for adults CT-11731

- tube- 1 piece
- breathing bag — 1 piece
- straight connector — 1 piece
- tube outlet and inlet valve with Luer-lock port — 1 piece
- tube 90 cm — 1 piece

Each design variant of the breathing circuit is supplied with instructions for use.

There are no **accessories** to the medical device. All components of the medical device are included in the package. If any part of the medical device is missing or damaged, the product should be replaced and discarded.

11. Technical characteristics

The composition of the breathing circuits includes (composition differs depending on the design variant):

1. Parallel connector
2. Straight connector 22M/15F
3. Y-shaped connector
4. Tube outlet and inlet valve with or without Luer-lock port
5. Heated wire
6. Tube of various diameters
7. Breathing latex or latex-free bag
8. Water collector
9. Automatic feed chamber
10. Hook
11. Straight connector 15F/22M
12. Connector 15F
13. Connector 22M/22F
14. Tube clamp
15. Heated inspiration or expiration tube

- 16. CO2 tube
- 17. Tube clamp
- 18. Nebulizer

Breathing tubes for circuits are produced in 3 diameters:

- For adults (ø 22 mm)
- For children (ø 15 mm)
- For newborns (ø 10 mm)

Technical characteristics:

Resistance to air flow

The air flow resistance without the breathing tube is 0.66 cm H₂O, when the breathing tube is connected to the connecting tube, the value is 0.94 cm H₂O.

The difference is 0.28 cm H₂O, which is equal to 0.028 kPa. The standard requirement for air flow resistance is 0.2 kPa. Thus, the Besmed breathing tube meets the standard requirements.

Safety of attachment of a connector with a normal size to a suitable size

The force used in the dynamometer is 61 N. The connector remains attached to the smooth end of the breathing tube, thus the device meets the requirements of this test.

Safety of connector connection to the breathing tube

The force used in the dynamometer is 66 N. The adapter is still attached to the breathing tube connector. Therefore, the product meets the requirements of this test.

Leakage test

Leakage from the breathing circuit should not exceed 150 ml/min (15.2 kPa·L/min) at an internal pressure of 30 hPa (30 cm H₂O).

The requirement of internal air pressure generation was fulfilled by an air generator. The pressure gauge controlled the air pressure of 6 kPa (equal to 45 mm Hg). The flow rate is 10 ml/min. The standard set by EN 12342 is 50 ml/min. Therefore, the breathing circuit has passed the leakage test.

Increase of the flow resistance with bending

The air flow resistance with an elongated breathing tube is 0.94 cm H₂O. When the bent breathing tube is connected to the connecting tube, the reading is 1.03 cm H₂O. The difference is 0.09 cm H₂O, which is equal to 0.009 kPa. This value is much less than 150% of the resistance value of the resistant breathing tube, so the device meets the requirements of this test.

Compliance verification

The fan provided the required internal air pressure. The pressure of the pressure gauge was controlled by an air pressure of 6 kPa (equal to 45 mm Hg). The volume of air required to maintain internal pressure is 10 ml. Tube compliance is expressed in milliliters per 1 kPa per meter of tube length. Breathing tube length is 1.2 meters. Therefore, this value is equal to 1.4 mL/kPa/meter

Breathing circuits and breathing mounts intended for use with flammable anaesthetic gases should meet the requirements of ISO 2878 during the ISO 2878 test.

The breathing circuit or its components marked as "antistatic" or "conductive" should have a maximum electrical resistance of 1.0 MOhm when tested in accordance with ISO 2878.

When evaluating the pipeline, the resistance limit should be 1 MOhm/m.

Breathing latex bag

- Size 0.5 L, 1 L, 2 L and 3 L
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

- Breathing bag connection hole, if provided, should be compatible with the breathing bag conforming to ISO 5362 and the breathing hoses (tubes) conforming to ISO 5367. This connection should be within 20° of the vertical axis. The air bag connection hole should not be located on the side of the patient's inhalation and exhalation valves.
- should comply with ISO 5362

Water collector

- Connector for connectors 15 M (diameter 15 mm) or 22 M (diameter 22 mm)
- Volume 30 ml
- Weight 65 g
- The material should be transparent or translucent
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

Automatic feed chamber

- Inlet / outlet size: 22M / 22M (diameter 22 mm)
- Maximum filling capacity is 150 ml.
- Compression volume of the empty chamber is 325 ml
- Compression volume of the filled chamber is 175 ml
- Diameter of the heated surface 12.2 cm
- Maximum flow 180L/min
- Resistance (cm of water column) 0.60 at a flow of 30 L/min
- Compliance 0.32 ml/cm H₂O.
- Weight 125 g
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

Nebulizer

- Optimal particle size (2-5 microns)
- Connector diameter 18 mm
- Size 20cc (volume 20 ml)
- Removable nebulizer cover with leak protection
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

Heated wire

- should comply with IEC 60601-1-2, IEC 60601-1-6, IEC 60601-1
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- Length 110, 180 or 200 mm
- Weight 45 g (at length 110 mm), 55 g (at length 180 mm) and 60 g (at length 200 mm)
- Hose diameter 22M

Connectors and tubes

Tube type (rifled, smooth-bore or stretchable) and the connection of the connector 15 F (15 mm) or 22F (22 mm) depends on the breathing circuit design variant. Connection of the connector should comply with ISO 5356-1.

The material of the connectors should be transparent or translucent and should be resistant to anaesthetic gases vapors. Connectors should be light, but strong enough, non-toxic, light-permeable with a compliance of less than 10%. Breathing tube connectors may be provided with protrusions, pads, or other means to facilitate attachment and detachment. These means should be well rounded. The cross-sectional area of the connector lumen should be not less than the lumen of the tube for installation in which it is intended. The minimum internal diameter of a curved shaped connector should be at least 80% of the nominal diameter, and its cross-sectional area should not be reduced by more than 10%.

If the connectors for the holes of the breathing circuit components can be disconnected by the operator without use of tools, they should be a 22 mm covered conical connector with a 15 mm coaxial covering conical connector conforming to ISO 5356-1, or a special connector incompatible with ISO 5356-1-compliant connectors and breathing tubes conforming to ISO 5367.

The surface of the connectors and tubes should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

The breathing tube should pass the ISO 5367 leakage test and not exceed 50 ml/min.

Connectors should be designed to have a small amount of dead space. The clearance should be free from irregularities.

A smooth internal channel should prevent condensation, which reduces the risk of bacterial growth.

The basic dimensions of the patient side of the connector should comply with EN ISO 5356-1.

The inspiration connector hole is a 22 mm covered conical connector that meets the requirements of EN ISO 5356-1.

The connector should be marked with an arrow indicating the direction of the gas flow. The axis of the joint should be located horizontally or within an angle of 50° to the horizontal plane.

The breathing tube should comply with EN ISO 5367.

Heated tubes should comply with IEC 60601-1-2, IEC 60601-1-6.

Connectors for other openings used for special purposes (e.g. pressure measurement, gas sample return) should not be compatible with ISO 5356-1, ISO 5356-2, or ISO 80369-7 connectors. Hole connectors should be equipped with means of reliable closing while they are not in use. Closure means should not get detached from the component.

Hooks and clamps should be clean, free from burrs, black particles, mechanical damage, foreign inclusions.

12. Description of materials contained in key functional elements / blocks / parts of the medical device, including materials in contact with the human body

All materials from which the components are made meet the specification and are allowed for use. The device and packaging do not contain any components of animal origin.

No.	The component part for which the material is used	Material, grade	Manufacturer
1	Breathing circuit tube for adults for children for newborns heated inspiration tube heated expiration tube CO2 tube	Polyvinylchloride EV08EFP02 EVA EV 302 ASIA polymer Corporation LDPE NA248	Green Field Chemtech co. Ltd. Formosa Plastics Corporation USI Company
2	Straight connector 15F/22M	LDPE (low density polyethylene) LDPE NA248	USI Company
3	Straight connector 22M/15F	Polypropylene 5090T	Formosa Plastics Corporation
4	Y-shaped connector	Polypropylene 5090 T	Formosa Plastics Corporation
5	Latex breathing bag	Latex 261419	laizhou huanhai latex technology co. ltd
6	Breathing bag	Polychloroprene CAS 25067-95-2 green pigment Green/980279	Nolato MediTor AB
7	Parallel adapter	Polypropylene 5090T	Formosa Plastics Corporation
8	Water collector	Polypropylene 5090T Polycarbonate S-2001R	Formosa Plastics Corporation Mitsubishi Chemical

			Holdings Corpora- tion.
9	Hook	Polypropylene 5090T	Formosa Plastics Corporation
10	Heated wire F&P	Polypropylene 5090T	Formosa Plastics Corporation
11	Parallel connector	Polypropylene 5090T	Formosa Plastics Corporation
	Automatic feed chamber	O-ring Silicone KE-1950-70- A The bottom plate of the chamber Aluminium AA1070 Float Silicone KE-1950-70- A Body Polycarbonate CP1- HT- PAL, Polypropylene 5090T Chamber cover Plastic LDPE NA248	Shin-Etsu Silicones Limit Taiwell Aluminum Corporation Shin-Etsu Silicones Limit Covestro (Hong Kong Limited Formosa Plastics Corporation USI Corporation
12	Nebulizer	Body upper part Plastic LDPE NA248 blue pigment Pigment No. Blue 15:0: CAS No, 147-14-8 Body lower part Polypropylene 5090T Tube Silicone KE-1950-70- A	USI Corporation Formosa Plastics Corporation Shin-Etsu Silicones Li
13	Connector 15F	Polypropylene 5090T	Formosa Plastics Corporation
14	Connector 15F/22M (inlet and outlet valve 15F/22M) without Luer-lock port with Luer-lock port	LDPE (low density polyethylene) LDPE NA248 blue pigment Pigment No. Blue 15:0: CAS No, 147-14-8	USI Company
15	Tube clamp	Polypropylene 5090T green pigment Pigment No. Green 7 : CAS No. 1328-53-6	Formosa Plasti Corporation
16		Polypropylene 5090T	

	Connector 22M/22F		Formosa Plastics Corporation
17	Package	HDPE polyethylene 8050	Formosa Plastics Corporation

13. Maintenance

Breathing circuits have a fully finished design and do not require any maintenance or repair interventions during operation.

14. Safe use conditions

-The procedure should be performed only by qualified medical personnel according to the instructions.

15. Order and preparation for work

15.1 Preparation for work

ATTENTION: Strictly follow aseptic technique and general precautions and procedures established by the protocols of the medical institution.

The procedure is performed according to the guidelines for medical professionals approved in a medical institution.

15.2 Operational procedure

15.2.1 Select the breathing circuit appropriate to the patient's age and clinical situation.

15.2.2 Open the package, remove the device, and check for damage. If the device is damaged, replace it with a new one.

15.2.3 When using a stretchable circuit, stretch it out to obtain necessary length by pulling the ends of the circuit.

15.2.4 Check the circuit for leaks.

15.2.5 Connect the circuit to the ventilator/patient's endotracheal tube (tracheostomy tube, laryngeal mask, etc.), make sure that the contour is leakproof.

15.2.6 To conduct gas analysis/capnography, connect the analyzer line to the corner connector port and make sure that the connection is tight.

15.2.7 When using an active humidifier, it is recommended to use a breathing circuit with a moisture collector. In this case, moisture collectors should be placed at the lowest point of the circuit branches. It is recommended to constantly monitor the liquid level, and, as necessary, empty the moisture collectors.

15.2.8 It is recommended to replace the device every 48 hours of ventilation or more often as contamination occurs.

15.2.9 The use of one respiratory circuit in several patients is not allowed.

15.2.10. Disinfection, sterilization and reuse are not provided, dispose of the device after use.

Note

- Device is used under the physician's supervision.
- If damage of the package is detected, the use of the breathing circuit is not recommended.
- Before using the breathing circuit, make sure that it is tight and that there is no leakage of the breathing mixture.
- Obstruction of the device by the discharge from the patient's airways, which disrupts ventilation and requires immediate replacement of the circuit, is possible.
- Fluid that accumulates in the breathing circuit can make it difficult to perform artificial lung ventilation. It is recommended to constantly monitor the liquid level, and, as necessary, empty the breathing circuit tubes and moisture collectors.
- It is not recommended to use an active humidifier or nebulizer in conjunction with a respiratory filter or with a respiratory filter with a heat-exchange device installed on the patient's connector, as this will lead to excessive accumulation of moisture in the circuit.

16. Cleaning and disinfection

Not provided. Single use device.

17. Environmental protection requirements, including requirements for safe destruction and disposal. Environmental protection.

- Requirements for the sanitary protection of air, soil and water bodies should correspond to the current sanitary and epidemiological rules and hygienic standards established in the country.

-The possibility of environmental pollution by production waste, generated in the process of production of the circuits, should be excluded.

- The technology of production of breathing circuits does not involve the formation of liquid waste.
- Disposal or destruction of technological losses resulting from the production of the circuits, as well as finished products, if it does not meet the requirements, should be carried out in accordance with current standards and requirements.
- When working with breathing circuits, no special measures to protect the environment from harmful effects are required.

Disposal:

- ☐ The used devices belong to class B waste according to SanPiN 2.1.7.2790.
- ☐ Systems that are not used for their intended purpose due to shelf-life expiration and other reasons are classified as class A, as epidemiologically safe waste and are disposed of as household waste.

18. Service life data

Finished breathing circuit sets meet the specifications and functional requirements for five years. Service life is not more than 48 hours.

19. Transportation and storage requirements

- ☐ Do not store under direct sunlight, at extreme temperatures, or high humidity.
- ☐ Store and transport at 0 ~ 40 ° C. The shelf life is 5 years.

20. Warranty

During the warranty period (within 5 years), the manufacturer can repair the device, if its failure is not caused by human actions.

21. Medical device labeling

21.1. Description of the MD labeling

Each consumer package should be provided with a printed product label, which should contain the following particulars:

- ☐ manufacturer's name and location (legal address, including country, and, in case it differs from the legal address, address(es) of production site(s));
- ☐ manufacturer's trademark (if available);
- ☐ device name (designation);
- ☐ package contents;
- serial number;
- ☐ inscriptions: "Non-toxic";
- ☐ inscription: "Do not use if the integrity of consumer packaging is violated";
- ☐ Lot number;
- ☐ warning about the need to check the integrity of the package before use;
- ☐ storage conditions;
- ☐ disposal conditions;
- ☐ registration certificate number;
- ☐ symbol "presence of instructions";
- ☐ symbol "about the need to read the instructions";
- symbol "single use";
- ☐ mark of conformity.

Note – The product label should be attached on the outside of the consumer package (packaged products should not come into contact with ink).

21.2. Description of the symbols, printed on the MD label

Batch code

Date of manufacture

Use before

Single use

Manufacturer

Authorized representative in the European Community

2460 Compliance with CE Directives


Do not use if the package is damaged


Device contains latex

Example of labeling for RF:

Breathing anaesthetic circuit in design variants

Breathing anaesthetic circuit for adults CT-11731

REF CT-11731

SN SO09040087

LOT SO09040087

2020/5/7

2025/5/7

Non-toxic

Besmed Health Business Corp.
No. 5, Lane 116, Wu-Kong 2nd Rd, Wu-Ku District,
New Taipei City, 248, Taiwan.
Tel: +886-2-2290-3959 | Fax: +886-2-2299-9076
www.besmed.com

PY M

2460


0°C 40°C


Check the package integrity before use

Breathing anaesthetic circuit in design variants

Breathing anaesthetic circuit for adults CT-16014

REF CT-16014

SN SO09040087

LOT SO09040087

2020/5/7

2025/5/7

Non-toxic

Besmed Health Business Corp.
No. 5, Lane 116, Wu-Kong 2nd Rd, Wu-Ku District,
New Taipei City, 248, Taiwan.
Tel: +886-2-2290-3959 | Fax: +886-2-2299-9076
www.besmed.com

PY M

2460

0°C 40°C

Check the package integrity before use

Labeling of secondary and transport packaging:

		Breathing anaesthetic circuit in design variants	
REF	CT-11731 Breathing anaesthetic circuit for adults CT-11731		
			
(01)04716770XXXXXX(10)MOXXXXXXX(17)XXXXXX			
LOT	MOXXXXXXX	 20XX-XX	 20XX-XX
 Besmed Health Business Corp. No. 5, Lane 116, Wu-Kong 2nd Rd, Wu-Ku District, New Taipei City 248, Taiwan			
Authorized representative of the manufacturer: Odyssey Medical LLC, Ap. 4, Build. 36, Zoya and Aleksander Kosmodemyansky st., Moscow, 125130			
       2460			

21.3. Sterilization and packaging

- The device(s) are packaged in a suitable rigid box. Package contents should be protected from impact. The packaging is in excellent condition: no breaks, holes or wet areas.
- The package should be securely closed. For this, a sufficiently strong polyethylene package is used.

The device is not sterile

There are no medicines in the composition of the medical device.

22. List of applicable national and international standards to ensure the safety, efficiency and quality of the medical device

Designation	Description
EN ISO14971:2012	Medical devices. – Application of risk management to medical devices
EN ISO 13485: 2012	Medical devices – Quality management systems. Requirements for regulatory purposes
EN 1041: 2008	Information supplied by the manufacturer of medical devices
EN 12342: 1998+ A1: 2009	Breathing tubes intended for use with anaesthetic apparatus and ventilators
ISO 5356-1:2015	Anaesthetic and respiratory equipment. Conical connectors. Part 1. Cones and sockets
ISO 80369-7:2016	Connecting parts with 6% taper (Luer) for syringes, needles and other medical equipment. Part 2. Luer tips
BS EN 60601-1-2:2015	Medical electrical equipment. Part 1-2. General requirements for basic safety and essential performance. Collateral Standard. Electromagnetic disturbances. Requirements and tests
IEC 60601-1-6:2010+AMD1:2013 CSV	Medical electrical equipment. Part 1-6. General requirements for basic safety and essential performance. Collateral standard. Usability
IEC 60601-1:2005/AMD1:2012/COR1:2014	Medical electrical equipment. Part 1.
DIN EN ISO 5367:2015-02	Anaesthetic and respiratory equipment. Breathing sets and connectors
ISO 5362:2006	Anaesthetic reservoir bags
DIN EN ISO 5356-2:2013-03	Anaesthetic and respiratory equipment. Conical connectors. Part 2. Screw-threaded weight-bearing connectors

Authorized representative of the manufacturer in the RF (name, address, telephone, mail)

Full or abbreviated company name (in accordance with the Constituent documents)	Odyssey Medical Limited liability company (Odyssey Medical LLC)
Legal address (in accordance with the Constituent documents)	Apt. 4, 36, Zoyi i Aleksandra Kosmodemyanskiki Moscow, 125130
Telephone (at the actual address)	8-925-145-61-34
E-mail	info@odyssey-medical.ru

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