
Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879


Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby swear (or affirm) that the attached documentation is true and unaltered copies of the original document.

Name: Linda Julia Davila
Title: Regulatory Affairs Specialist

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
Pulmonary Artery Vent Cannula

Product Description

The cannula consists of a flexible plastic tube that terminates in a 0.64 cm (1/4 in) ribbed female luer. A soft, pliable, tapered tip with 12 side ports and a 12 ga introducer needle aid in the insertion of the cannula. A moveable collar acts as a depth marker.

Sterile, nonpyrogenic, single use.

Indications for Use

This cannula is intended for use in venting the pulmonary artery during cardiopulmonary bypass surgery for up to 6 hours or less.

The pulmonary artery vent cannula (PAVC) is a simple, effective, and atraumatic technique for decompression of the left side of the heart during aortic valve replacement, mitral valve repair or replacement, or coronary artery bypass.

The PAVC decreases the blood flow to the left side of the heart from bronchial flow, pulmonary vasculature, and spillover from the right side of the heart.

The PAVC also minimizes the possibility of air in the left heart during coronary artery bypass surgery.

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Warnings

Gravity drainage is suggested. In procedures where relatively high flows may be encountered, it is important to monitor the degree of suction to reduce the possibility of air embolism.

Precautions

The catheter is made of PVC, which will stiffen at hypothermic cardioplegia temperatures. Necessary care must be taken when manipulating the heart.

Note: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, the type of surgical procedure, and the type of delivery and suction systems utilized.

When removing the cannula, withdraw it completely before tightening the purse string suture to prevent catching the suture in the perforated section of the cannula.

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Applicable where phthalates are identified on product label: Due to the presence of phthalates in the product, the clinician must weigh the medical benefits of product use against the drawbacks of phthalate exposure for male children and pregnant or nursing women.

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Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique. Remove the flexible tip protector and discard it.
2. Withdraw the introducer needle from the cannula completely to insure that the introducer needle can be withdrawn smoothly. Inspect the point on the introducer needle to insure that the point has not been inadvertently bent or damaged. Carefully replace the introducer needle in the cannula and turn the locking fitting to secure the introducer needle in place. Do not overtighten.

Note: Use caution when reinserting introducer needle.

3. Place a purse string suture at the desired cannulation site on the pulmonary artery (with or without pledgets). Draw the suture ends through a tourniquet.
4. Insert the needle into the pulmonary artery through the center of the suture. The cannula will act as its own dilator. Rotate and insert the cannula until all of the side holes are positioned in the pulmonary artery (Figure 2).

5. Loosen the introducer needle by rotating the locking fitting counterclockwise. Withdraw the introducer needle and continue to insert the cannula to the depth collar. Clamp or pinch the cannula at a convenient point on the flexible tubing below the inlet fitting using your fingers or a suitable tubing clamp (not supplied).
6. Secure the tourniquet and tie it to the cannula at a convenient site above the collar.
7. Connect the cannula to a suction line from the extracorporeal circuit.
8. After the procedure, withdraw the cannula completely and secure the purse string suture. Dispose of the cannula according to hospital policy.

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Disclaimer of Warranty

ALTHOUGH THE PULMONARY ARTERY VENT CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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Disclaimer of Warranty

ALTHOUGH THE PULMONARY ARTERY VENT CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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Medtronic

Europe

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl

Route du Molliou 31

Casa Postale 84

CH - 1131 Tolochenaz

Switzerland

Internet: www.medtronic.co.uk

Tel. 41-21-802-7000

Fax 41-21-802-7900

Authorized Representative in the European

Community

Medtronic B.V.

Eari Bakkenstraat 10

6422 PJ Heerlen

The Netherlands

Tel. 31-45-566-8000

Fax 31-45-566-8668

Asia-Pacific

Japan

Medtronic Japan

Comodio Shiodome 5F

2-14-1 Higashi-Shimbashi, Minato-ku

Tokyo 105-0021

Japan

Tel. 81-3-6430-2011

Fax 81-3-6430-7140

Australia

Medtronic Australasia Pty. Ltd.

Business Address:

97 Waterloo Road

North Ryde NSW 2113

Mailing Address:

PO Box 945

North Ryde NSW 1670

Australia

Tel. 61-2-9857-9000

Fax 61-2-9878-5100

Asia

Medtronic International Ltd.

Suite 1602 16/F, Manulife Plaza

The Lee Gardens, 33 Hysan Avenue

Causeway Bay

Hong Kong

Tel. 852-2891-4068

Fax 852-2591-0313

Americas

Latin America

Medtronic Latin America

3750 NW 87th Avenue

Suite 700

Miami, FL 33178

USA

Tel. 305-500-9328

Fax 786-709-4244

Canada

Medtronic of Canada Ltd.

5733 Kilmat Road

Mississauga, Ontario L5N 1W3

Canada

Tel. 905-826-6020

Fax 905-826-6820

Toll-free in Canada:

1-800-268-5346

United States

Manufacturer:

Medtronic, Inc.

710 Medtronic Parkway

Minneapolis, MN 55432

USA

Internet: www.medtronic.com

Tel. 763-514-4000

Fax 763-391-8100

Toll-free in the USA:

1-800-328-2518

(24-hour consultation service)

M940825A001


© 2003, 2010 Medtronic
M940825A001 Rev. A

Medtronic

DLP®

Left Heart Vent Catheter
Cathéter de décompression du ventricule gauche
Linksherz-Entlüftungskatheter
Catéter para ventilación del lado izquierdo del corazón
Ontluchtingscatheter voor de linkerhart helft
Catetere ventricolare sinistro
Ventileringskateter, venstre side av hjertet
Sydämen vasemman puoliskon ilmanpoistokatetri
Kateter för avluftning av vänsterhjärtat
Καθετήρας εξάερωσης αριστερής κοιλίας
Udluftningskateter til venstre hjertekammer
Cewnik do odpowietrzania lewego serca
Cateter de ventilação do coração esquerdo
Sol Kalp Vent Kateteri
Odvzdušňovací katétr pro levé srdce
Bal kamrai szívókatéter
Дренажный катетер для левых отделов сердца

Handwritten signature
May 29, 2018

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
Bruksanvisning • Οδηγίες χρήσης • Brugsanvisning • Instrukcja użytkowania
Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati utasítás
Инструкция по эксплуатации

(USA) Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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Left Heart Vent Catheter

Product Description

The Silicone Left Heart Vent consists of a flexible silicone tube with a perforated 8.9-cm (3.5-in) distal segment. A small tapered silicone tip is attached to the perforated distal end to aid in the insertion of the catheter across the mitral valve. A malleable stainless steel introducer is provided with the catheter. Depth markings on the tube indicate insertion depth.

The PVC Left Heart Vent consists of a flexible plastic tube with a perforated distal segment. A small tapered silicone tip is attached to the perforated end to aid in the insertion of the proximal end of the catheter across the mitral valve. Depth markings on the tube indicate insertion depth.

The Malleable PVC Left Heart Vent consists of a flexible plastic tube with a perforated distal segment. A malleable thin wire within the catheter tube permits precise shaping. Depth markings on the tube indicate insertion depth.

The PVC Left Heart Vent Catheter with Pressure Monitoring Line consists of a one-piece, dip-molded body with side holes in the tip and terminates in a 0.95-cm (3/8-in) nonvented connector. The proximal end features a pressure monitoring line with a 3-way stopcock. Depth markings beyond the last set of holes in the tip aid in positioning. Overall length of the catheter is approximately 43.8 cm (17.25 in).

Sterile, nonpyrogenic, single use.

Indications for Use

This catheter is intended for use in venting the left heart during cardiopulmonary bypass surgery up to six hours or less.

Contraindications

This device is not intended for use except as indicated above.

Warnings

Gravity drainage is suggested in procedures where the left heart is kept closed in order to minimize the possibility of air embolism and over-distention of the left ventricle.

Gravity drainage is not suggested where relatively high flows may be encountered.

Precautions

Note: Proper surgical procedures and techniques are the responsibility of the medical profession. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, the type of surgical procedure, and the type of suction systems utilized.

The PVC models will stiffen at cold temperatures. Necessary care must be taken when manipulating the heart to avoid inadvertent trauma to delicate cardiac tissues.

When removing the catheter, withdraw it completely before tightening the purse string suture to prevent catching the suture in the perforated distal tip of the catheter.

Please tighten luer caps before using.

The catheter is designed and intended for single use only. **DO NOT REUSE.**

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events, and dislocation. Vessel damage and complications at the puncture site may occur if the instructions for Use are not followed.

Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the catheter onto the sterile field utilizing an aseptic technique.
2. Prior to insertion, bend the catheter into the desired shape. If the catheter has an introducer, simply form the catheter with the introducer in place. If using a PVC model without an introducer, bend the catheter into the desired shape and hold in cold sterile saline to temporarily form it to the new shape. If using the Malleable PVC model, make certain that the wire within the catheter body is on the inside of the formed angle to avoid kinking or occluding the lumen of the catheter. Avoid forming abrupt angles.
3. Retract the right atrium and place a purse string suture at the cannulation site in the right superior pulmonary vein adjacent to its junction with the left atrium. Draw the suture through a tourniquet/ligature tube. Make an incision within the purse string suture.

PVC Left Heart Vent Catheter with Pressure Monitoring Line models: Prime the pressure monitoring lumen of the catheter with heparinized saline and close the stopcock toward the lumen to retain the heparinized saline in the lumen.

4. Moisten the catheter in sterile solution before inserting. Insert the tip of the catheter just through the wall of the right superior pulmonary vein and advance into the left atrium.

Note: It is recommended that the aorta be cross clamped prior to the insertion of the vent catheter.

5. Insert the catheter at least past the last set of holes on the distal tip of the catheter or at least to the first line on the catheter. During insertion, palpate the left ventricle to insure proper positioning of the catheter across the mitral valve.
6. Tighten the tourniquet and secure to the vent catheter at a convenient site. Remove the introducer (if included). Connect the vent catheter to the suction tubing provided in the perfusion circuit.

PVC Left Heart Vent Catheter with Pressure Monitoring Line models: Attach a heparinized, saline-filled syringe to the female luer that is perpendicular to the pressure monitoring line tubing axis. Aspirate the pressure monitoring line until blood appears in the lumen. Tap the syringe to force any air drawn up during aspiration to the top of the syringe. Inject 2-3 cc of heparinized saline to fill the pressure monitoring line.

Attach a preprimed pressure monitoring line from the pressure monitoring device to the 3-way stopcock. Make certain the stopcock is open to the pressure monitoring line of the catheter. Balance and zero the transducer.

7. After the procedure, remove the vent catheter from the heart and tie the purse string suture (see Precautions).

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE LEFT HEART VENT CATHETER, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE LEFT HEART VENT CATHETER, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Technical characteristics:

1. Cannula length: 15.2 cm
2. Effective length of the cannula: 12.3 cm
3. Outer diameter of the flexible tip is 16 Fr (5.3 mm)
4. Diameter of the flexible tube:
Outer: 5.3 mm
Inner: 4 mm
5. Length of protection of the flexible tip: 4.5 cm
6. Needle gauge is 12G (2.7mm)
7. Needle length: 17.8 cm
8. Outer diameter of the needle, max.: 2.8 mm
9. Diameter of the needle limiter: 3.3 mm
10. Length of the needle limiter: 9 mm
11. Number of lateral ports: 12 pieces
12. Vent slot
The maximum diameter is 2 mm
13. Product total mass: 6.28 g
14. Diameter (Locking male luer fitting), mm:
Outer 8.0 ±0.1

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
For Left Heart Vent Catheters DLP:

an example of the individual package marking of the model 12004

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
Technical characteristics:

1. Catheter length: 38.1 cm
2. Catheter effective length: 100 mm
3. Diameter of the flexible tube:
Outer: 3.3 mm($\pm 0,2$ mm)
Inner: 2.3 mm
4. Slip connector diameter: 6.4 mm
5. Vent slot diameter: 1.7 mm
6. Number of vent slots: 9 pieces
7. Vent port diameter:
Inner (Locking female luer fitting): 4.270-4.315 mm
Outer (Locking female luer fitting): 5.7 mm

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

Technical characteristics:


1. Catheter length: 38.1 cm
2. Catheter effective length: 100 mm
3. Diameter of the body:
Outer: 4.3 mm ($\pm 0,2$ mm)
Inner: 2.8 mm
4. Slip connector diameter: 6.4 mm
5. Vent slot diameter: 2.4 mm
6. Number of vent slots: 9 pieces
7. Vent port diameter:
Inner (Locking female luer fitting): 4.270-4.315 mm
Outer (Locking female luer fitting): 5.7 mm
8. Product total mass: 6.23 g
9. Minimum force at break is 15 N
10. Distal Perforation Length: 10,4 mm
11. Distal End to Depth Marks Lengths: (20,3 – 30,5 – 40,6 – 50,8 – 61,0 – 71,1 – 81,3 – 91,4 – 101,6) mm
12. Radius of Rounding of Distal Tip: 2,16 mm

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

Technical characteristics:

1. Catheter length: 32 cm
2. Stylet length: 34.5 cm
3. Catheter effective length: 100 mm
4. Outer diameter of the tip: 13 Fr (4.3 mm)($\pm$ 0,2 mm)
5. Diameter of the flexible tube:

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Technical characteristics:

1. Catheter length: 33 cm
2. Catheter effective length: 100 mm
3. Outer diameter of the tip: 13 Fr (4.3 mm) ($\pm 0,2$ mm)
4. Diameter of the flexible tube:
Outer: 4.3 mm ($\pm 0,2$ mm)
Inner: 2 mm
5. Slip connector diameter: 6.4 mm
6. Vent slot diameter: 2 mm
7. Number of vent slots: 14 pieces
8. Vent port diameter:
Inner (Locking female luer fitting): 4.270-4.315 mm
Outer (Locking female luer fitting): 5.7 mm
9. Product total mass: 11.74 g
10. Minimum force at break is 15 N
11. Distal Perforation Length: 21, 6 mm

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Description:

Transport package is made of corrugated cardboard. The sizes of a transport package vary depending on the order.

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Scope of delivery:**For left heart vent catheters DLP, variants:**

The individual package contains (for models (12116, 12118, 12110, 12113, 12115, 12001, 12008):

-Catheter

The individual package contains (for models ((12101, 12016, 12002, 12220):

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2. Application instruction (with a instruction leaflet) - 1 piece;
The number of individual boxes in the transport package varies depending on the order.

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Information in Russian sticker

The following information should be applied on a Russian-language sticker (word order and text font may vary)

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following information:

Name of medical device in Russian language;
Name of the manufacturer;
Country of origin;
Code UPN (unique product number);
Code CFN (customer facing number);
Date and number of Registration Certificate;
Sign of conformity in Russia;
Warning and informing symbols and inscriptions (if applicable).
Exsample of the sticker

	Name of medical device	
	Manufacturer: Medtronic Inc.	Code UPN
	Country of origin: USA	Code CFN

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All Rights Reserved