Handwritten signature

Henke-Sass, Wolf GmbH
Keltenstrasse 1 · 78532 Tuttlingen · Germany

Henke-Sass, Wolf GmbH
Keltenstrasse 1
78532 Tuttlingen · Germany
Telefon: +49 7462 9466-0
Telefax: +49 7462 9466-5000
info@henkesasswolf.de
www.henkesasswolf.de

Certification for Information in Dossier

The submitted dossier for the Russian application contains the following:

- Addendum of Instructions for use of **Rigid medical endoscopes**

I certify that the documents which are submitted are true and accurate to the best of my knowledge. In addition, Henke-Sass, Wolf GmbH takes responsibility for any legal obligations related to this submission.

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EN

Table of contents

1.	About this document	12
1.1	Purpose	12
1.2	Symbols used in this document	12
2.	Intended use	13
2.1	Intended use – sheaths	13
2.2	Intended use – adaptors	13
	Intended User	13
2.3	Contraindication	13
3.	Safety information	13
4.	Guideline conformity	14
5.	Description	14
5.1	Markings	14
5.2	Available models and combinable products	14
6.	Preparation for use	15
6.1	Visual inspection and function check	15
6.2	Provisioning	15
7.	Use	15
8.	Processing	16
8.1	Warning and precautions	16
	Limitations on reprocessing	16
8.2	Cleaning and disinfection	16
	Transport	16
	Manual cleaning / pre-cleaning and chemical disinfection	16
	Disinfection	17
8.3	Sterilisation	17
	Fractionated pre-vacuum method ¹⁾	18
	Packaging	18
8.4	Tray configuration	18
9.	Assembly	20
9.1	Sheaths and obturators	20
9.2	Adaptor with stopcock	20
10.	Disassembly	20
11.	Storage	20
12.	Ambient conditions	21
13.	Service and maintenance	21
14.	Disposal	21

1. About this document

1.1 Purpose

This document describes the correct handling and function of the product, as well as recommended processing methods.

This document may not be used to carry out examinations or surgeries, nor may it be used for training purposes.

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If you as the user of this product believe that you require more detailed information regarding the product's use and maintenance, please contact your representative.

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	Instructions for preventing personal injury
	Instructions for preventing material damage
	Information to facilitate understanding or workflow optimisation
	Prerequisite
	Instruction

2. Intended use

Zimmer Biomet has a range of arthroscopic instruments used during surgical and diagnostic procedures. The arthroscopic instruments include sheaths, obturators, trocars and adaptors.

- Carefully read these instructions before using the products. Keep them in a safe place for future reference.

2.1 Intended use – sheaths

These products are intended for the arthroscopic diagnosis and treatment of the skeletal system.

2.2 Intended use – adaptors

The adaptors have the same intended use as the sheaths, as they are always used in combination.

Intended User

The products must be used only in medical facilities by trained and skilled medical personnel. The products must not be used if according to a qualified physician, the general condition of the patient is not adequate or if the arthroscopic methods are contraindicated.

2.3 Contraindication

- Do not use the device if it has been already used to treat patients with suspected or verified Bovine Spongiform Encephalopathy, Creutzfeldt-Jakob disease / Variant Creutzfeldt-Jakob disease.

DANGER

Patient hazard!

- Treat surgical patients identified as at-risk for Creutzfeldt-Jakob disease (CJD) and related infections with single-use instruments.
- Dispose of products that have been used or are suspected of having been used on a patient with CJD after surgery according to current national recommendations.

3. Safety information

The instrument may only be used by trained medical professionals in medical facilities.

- When the instrument is delivered, inspect the delivery for completeness and damage.
- Read, observe and keep the instructions for use.
- Use the instrument only as intended, see "Intended use" on page 13.
- When endoscopic instruments are used together with other electromedical devices, the total patient current leakage may be additive. Ensure BF (CF) conditions are observed when using the endoscopic instruments with electromedical devices.
- Failure to use protective filters or suitable filtering spectacles during the activation of a surgical laser beam may cause eye damage to the user.
- The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the member state in which the user and/or patient is established.
- Ensure that any modifications or repairs are carried out by persons authorised by Henke-Sass, Wolf. No modifications of this equipment are allowed by the user, or the warranty is void.

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- Ensure that the processing, material and personnel are suitable for achieving the results necessary.

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The instruments are delivered non-sterile as reusable products.

- Clean / disinfect and sterilise the instrument prior to initial use as well as each subsequent use of the instrument.
- Bring the instrument to the decontamination area after use. Observe valid protective measures to prevent contaminating the environment.

**WARNING**

Risk of injury due to a faulty instrument!

- Carry out visual inspection and function check prior to each use.
- Only use instruments which are in perfect condition.
- Observe valid local operator regulations for all manual cleaning and drying processes.

- During storage, transport and processing, ensure that the instrument is not subjected to mechanical strain.

- Observe the following measures to prevent material damage:
 - Do not subject the instrument to impact.
 - Do not bend the sheath.
 - Put the instrument down carefully.
 - Transport instruments individually and store them safely by using a screen basket or container.

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5. Description

5.1 Markings

- Article number
- LOT number
- CE mark with identification number of the notified body where applicable:
Instrument conforms to the requirements of the guideline 93/42/EWG
- Writing
 - ZIMMER BIOMET
 - HSW

5.2 Available models and combinable products

Information about available models and combinable products can be requested at Zimmer Biomet.

6. Preparation for use

6.1 Visual inspection and function check

WARNING

Risk of injury due to a faulty instrument!

- Carry out a visual inspection and function check prior to initial use as well as each additional use.
- Only use instruments which are in perfect condition.
- Clean / disinfect and sterilise the instrument prior to initial use as well as each subsequent use of the instrument.
- Ensure that no parts are missing or loose.
- Ensure that there are no residual cleaning agents or disinfectants on the instrument.
- Inspect the entire instrument, particularly the sheath, for contaminants and damage of any type, such as dents, scratches, cracks, bending and sharp edges.
- For instruments with locking device: inspect the locking device for contaminants and damage to ensure a fixed and secure connection.
- For instruments with stopcocks: inspect all components of the stopcock for function and damage.
- Inspect free passage of the working channels.
- Prepare the instrument for processing immediately after use to prevent surface drying of the contaminants.

6.2 Provisioning

- Clean / disinfect and sterilise the instrument prior to initial use as well as each subsequent use of the instrument, see "Processing" on page 16.

7. Use

For instruments with working channel:

- If an instrument is inserted (if there is a stopcock on the working channel): do not move the lever on the stopcock because otherwise the instrument could be damaged or sheared off.
- Prepare the instrument for processing immediately after use to prevent surface drying of the contaminants.

EN

8. Processing

8.1 Warning and precautions

- Observe the country-specific regulations and laws concerning cleaning of instruments for patients with Creutzfeldt-Jakob-Disease, CJK-on-spec or its possible variants, Bovine Spongiform Encephalopathy or Transmissible Spongiform Encephalopathy.

Limitations on reprocessing

- Using multiple sterilization methods may significantly reduce the performance of the instrument and is not recommended.
- Proper processing has a minimal effect on this instrument. Lifetime is normally determined by wear and damage due to use.
- Damage incurred by improper processing will not be covered by the warranty.

8.2 Cleaning and disinfection

Transport

- Always store instrument securely and transport it to processing in a closed container to prevent damage to the instrument and the contamination of the environment.

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- Do not use fixating cleaning agents or hot water ($> 40^{\circ}\text{C}$) as it can cause fixation of the contaminants and jeopardise successful cleaning.
- Do not scratch contaminants off with hard objects as this can cause damage to the instrument.

- For instruments with stopcocks: All components of the stopcocks are dismantled, see chapters 9.2 and 10.
- Remove coarse contamination from the instrument. Clean the instrument with a soft brush under cold tap water until all visible contaminants have been removed.
- For working channels, holes and threads: Rinse with a water pistol for at least 10 seconds at a pressure of at least 3.8 bar (absolute).
- Disinfect the instrument. In doing so, observe the specifications of the disinfectant solution manufacturer regarding temperature, concentration and application time.

Non-compliance with the manufacturer's specifications can result in damage to the instrument.

- Rinse instrument with deionised water.
- Dry working channel with compressed air.
- Dry instrument with lint free cloths.
- Carry out visual inspection, function check and servicing.

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The instruments are suitable for prevalent machine methods of cleaning and thermal disinfection. In doing so use gentle cycles for instruments and suitable cleaning agents and disinfectants. The instructions of the machine, cleaning agent and disinfectant manufacturers must be observed. The cleaning and disinfectant result must be confirmed by the machine, cleaning agent and disinfectant manufacturers in cooperation with the user.

- Fix the instruments to the corresponding tray in such a way that a successful cleaning procedure is possible and that damage is prevented during cleaning (see pictures of the trays in chapter 8.4).
- Place the tray in an automated washer.
(The automated cleaning was validated with the cleaning unit Miele Professional G 7836 CD. For use with another washer/disinfector make sure that it is also in compliance with EN ISO 15883-2.)

- Program the washer with the following parameters and activate the wash.

Step	Process step	Reagents	Time	Temperature
1	Pre-cleaning	Tap water	3 min	Cold tap water ($\leq 23\text{ }^{\circ}\text{C}$)
2	Drain	–	–	–
3	Cleaning	Tap water, • Dosing: 0.5% neodisher MediClean forte at $55\text{ }^{\circ}\text{C}$	5 min	$55\text{ }^{\circ}\text{C}$ (set point)
4	Drain	–	–	–
5	Rinsing	Deionised water	1 min	Cold deionised water ($\leq 23\text{ }^{\circ}\text{C}$)
6	Drain	–	–	–
7	Rinsing	Deionised water	1 min	Cold deionised water ($\leq 23\text{ }^{\circ}\text{C}$)

Disinfection

Thermal disinfection has been validated using the following parameters.

- Time: 95 seconds
- Temperature: $95\text{ }^{\circ}\text{C}$
- Ensure that the exteriors of the instrument are dry. If necessary, dry with lintfree cloths.
- Dry working channels with compressed air.
- Carry out visual inspection, function check and servicing.

8.3 Sterilisation

Prior to each sterilisation, the instruments must be cleaned and disinfected according to the methods in these instructions for use.

- For instruments with stopcocks. Assemble all stopcocks prior to the sterilisation (see chapter 9.2).
- Double wrap the tray prior to sterilisation according to AAMI ST 79. Sterilise instruments in suitable packaging to prevent subsequent contamination.
(The validation was performed with Broemeda Amcro Flexibles sterilisation wraps $900 \times 900\text{ mm}$ (Broemeda item number 68170990).)

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Comply with specified process parameters. The specified parameters have been validated to ensure the sterility of the instruments. Deviating process parameters could damage the instrument. In this case, the guarantee and warranty shall become void.

- Sterilise instruments.
- When the sterilisation process has ended, allow the instruments to cool gradually to room temperature.

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¹⁾ The sterilisation validation was validated with the steam steriliser Selectomat HP 686-1HR of MMM.

²⁾ 18 minutes is the maximum exposure time that should be used for the prevacuum cycle.

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
➤ Pack the products into suitable sterilisation packaging systems (e.g. STERICLIN pouch used during sterilisation validation) acc. to ISO 11607 and/or AAMI / ANSI ST77 2006 in order to be sterilised.

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


9. Assembly

9.1 Sheaths and obturators

Arthroscopic sheaths and corresponding obturators are colour coded blue. To introduce an obturator into the sheath with the Quick-Lock-Device:

- Insert the obturator, lining-up the arrow markings on the sheath and the obturator.
- Lock in by gently pushing until the locking device clicks in.
- To remove the obturator, push the button to activate the click mechanism and release the obturator. Gently pull it out of the sheath.

- Never use force to insert obturator. It should glide in easily. Otherwise check for correct size or that it is correctly aligned.

9.2 Adaptor with stopcock

i

In order to ensure sterility, only use grease which is suitable for medical instruments for the stopcock (1) and stopcock plug (3 or 4).

- Apply a thin layer of grease to the contact surfaces of the stopcock plug (3 or 4).
- Mount the stopcock plug (3 or 4) and fix with the stopcock nut (2).
- Remove excess grease.
- Dismantle the stopcock for sterilisation.
- All the individual parts of the stopcock (1, 2, 3 or 4) are to be packaged for sterilisation in suitable sterilisation packaging systems (for example STERICLIN as at the validation of sterilisation) in accordance with ISO 11607 and/or AAMI /ANSI ST77 2006.

Adaptors can be attached in the same way as the obturators (quick-lock) with the alignment arrows.

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As long as instruments are stored unsterile in the original packaging, the following storage conditions apply

- Avoid direct sunlight.
- Store instrument either in the original packaging or in a screen tray/container.
- Ensure that the instrument is stored securely.
- Apply the respective valid national provisions when storing in a sterile condition.

12. Ambient conditions

	Operation	Transport and storage
Temperature	10 °C - 40 °C	-40 °C - 70 °C
Relative humidity	5 % - 75 %	5 % - 95 %
Atmospheric pressure	70 kPa - 106 kPa	70 kPa - 106 kPa

13. Service and maintenance

Henke-Sass, Wolf does not supply original parts to independent workshops or other manufacturers. Thus only Henke-Sass, Wolf is in a position to carry out repairs using original parts. The original technical specifications and the operational safety of the instrument can only be guaranteed by using original parts. The warranty for Henke-Sass, Wolf products shall become void if repairs are carried out by a workshop not authorised by Henke-Sass, Wolf. In this case Henke-Sass, Wolf is also no longer responsible for the technical specifications or safety of the product.

- Have the products repaired by Henke-Sass, Wolf only. For service, send the defective product to the address of the service centre. Please ask your representative for the address.
- Clean, disinfect and sterilise the product thoroughly prior to returning it for repair.
- Ideally, send in the product in its original packaging. If this is not possible, package the product carefully to secure it for transport.
- Henke-Sass, Wolf is not liable for damage resulting from improper shipping.

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EN

Table of contents

1.	About this document	14
1.1	Purpose	14
1.2	Symbols used	14
2.	Intended use	15
3.	Safety information	15
4.	Testing, handling and maintenance	16
5.	Guideline conformity	18
6.	Description	17
6.1	Construction	17
6.2	Markings	17
6.3	Available designs and sizes	17
6.4	Combinaible products	17
7.	Preparation for use	18
7.1	Visual inspection and function check	18
7.2	Provisioning	18
8.	Use	19
9.	Processing	19
9.1	Warnings and precautions	19
9.2	Cleaning and disinfection	19
	Manual cleaning / pre-cleaning and chemical disinfection	19
	Machine cleaning and thermal disinfection	20
9.3	Sterilisation	20
	Steam sterilisation (autoclaving)	21
	Fractionated pre-vacuum method ¹⁾	21
9.4	Special precautions: Pathogens of Transmissible Spongiform Encephalopathy	21
9.5	Processing restrictions	21
9.6	Tray configuration	22
10.	Assembly	23
11.	Disassembly	24
12.	Storage	24
13.	Ambient conditions	25
14.	Service and maintenance	25
15.	Disposal	25

1. About this document

1.1 Purpose

This document describes the correct handling and function of the rigid endoscope, as well as recommended processing methods.

This document may not be used to carry out endoscopic examinations or surgeries, nor may it be used for training purposes.

The respective current version of this document can be requested from Zimmer Biomet.

If you as the user of this endoscope believe that you require more detailed information regarding the product's use and maintenance, please contact your representative.

1.2 Symbols used

The following symbols are used to make it easier for you to access the information.

	Instructions for preventing personal injury
	Instructions for preventing material damage
	Information to facilitate understanding or workflow optimisation
	Prerequisites
	Instruction

2. Intended use

Rigid medical endoscopes are used to visualise body cavities. Each endoscope was developed for diagnostic and surgical procedures in one of the following fields of application:

- Arthroscope: arthroscopic procedures
- Sinuscope: sinusoscopic procedures
- Otoscope: otoscopic procedures
- Laryngoscope: laryngoscopic procedures
- Laparoscope: laparoscopic procedures
- Hysteroscope: hysteroscopic procedures
- Cystoscope: cystoscopic procedures

For the benefit and safety of patients, physicians must select a method which they consider suitable based on their experience. If you as the user of this endoscope believe that you require more detailed information regarding the product's use and maintenance, please contact your representative.

3. Safety information

The endoscope may only be used by trained medical professionals in medical facilities.

- After delivery, inspect the endoscope for completeness and damage.
- Read, observe and keep the instructions for use.
- Use the endoscope only as intended, see "Intended use" on page 15.
- When endoscopic instruments are used together with other electromedical devices, the total patient current leakage may be additive. Ensure BF (CF) conditions are observed when using the endoscopic instrument with electromedical devices.
- Failure to use protective filters or suitable filtering spectacles during the activation of a surgical laser beam may cause eye damage to the user.
- The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the member state in which the user and/or patient is established.
- Ensure that any modifications or repairs are carried out by persons authorised by Henke-Sass, Wolf. No modifications of this equipment are allowed by the user, or the warranty is void.

For storage, transport and processing, ensure that the endoscope is not subjected to mechanical strain, particularly to prevent damage to the sensitive lens system.

WARNING

Risk of infection to the patient or medical professionals!

The endoscopes are delivered non-sterile as reusable products.

The state of the art and national laws require the observance of validated processes.

In general, users are responsible for the validation of their processes.

- Ensure that the processing, material and personnel are suitable for achieving the results necessary.
- Observe any valid local operator regulations for all manual cleaning and drying processes.
- Clean / disinfect and sterilise the endoscope prior to initial use as well as each subsequent use of the endoscope.
- Bring the endoscope to the decontamination area after use. Observe valid protective measures to prevent contaminating the environment.

EN

**WARNING****Risk of burns!**

The optical fibres emit high-energy light at the distal end of the endoscope. This can cause the temperature of the body tissue to rise to 41 °C.

- Avoid direct contact of the distal end with body tissue or flammable materials as it can cause burns.
- Reduce the light intensity of the cold light source when working near body tissue or flammable materials.

**WARNING****Risk of injury due to faulty endoscopes!**

- Carry out visual inspection and function check prior to each use.
- Only use endoscopes which are in perfect condition.

4. Testing, handling and maintenance

Endoscopes are precision medical instruments, and handling them requires great care.

- Inspect the endoscope for damage prior to and after use.
- If the endoscope is damaged, discontinue use and contact the manufacturer.
- Do not subject the endoscope to impact.
- Put the endoscope down carefully.
- Hold endoscope only by the ocular funnel / main part and not by the sheath.
- Do not bend the sheath.
- Do not bend the body after inserting the endoscope into the body. A piece broken off the endoscope can become lodged in the soft tissue or no longer appear in the endoscope's field of vision and thus remain in the body.
- Transport endoscopes individually and store them safely by using a screen basket or container.

5. Guideline conformity


CE

CE 2797

The CE marking of the medical product complies with the guideline 93/42/EWG. It applies only when the product and/or packaging features this marking. In addition, products of Class IIa and higher are labelled with the identification number of the respective notified body.

6. Description

6.1 Construction

- 1 Distal end
- 2 Sheath
- 3 Main part
- 4 Ocular funnel
- 5 Proximal end
- 6 Irradiation surface of the illumination fibres
- 7 Connection for illumination fibre, type ACMI
- 8 Adapter for illumination fibre, type Wolf, pre-assembled
- 9 Adapter for illumination fibre, type Storz / Olympus (assembly, see "Assembly" on page 23.)
- 10 For endoscopes with working channel: Adapter with stopcock (disassembly, see "Disassembly" on page 24.)

6.2 Markings

- Article number
- Serial number
- CE mark with identification number of the notified body where applicable. Endoscope conforms to the requirements of the guideline 93/42/EWG
- For autoclavable endoscopes: Writing "autoclavable"
- Specification of the direction of view
- Writing GERMANY
- Writing  HWS
- For large-screen lenses: Writing HM
- For high-resolution optics developed specifically for the latest Full HD/4K camera generations: Writing HD4K

6.3 Available designs and sizes

The endoscopes are available in the following designs and sizes:

- Straight endoscopes
- Angled endoscopes
- Endoscopes with working channel
- Sheath diameter 1.9–11 mm
- Working channel diameter 5.1–6.4 mm

6.4 Combinable products

Endoscopes are generally compatible to camera systems and light sources that use interfaces according to ISO TS 18339. Information about available models and combinable products can be requested at Zimmer Biomet.

EN

7. Preparation for use

7.1 Visual inspection and function check

WARNING

Risk of injury due to faulty endoscopes!

- Carry out a visual inspection and function check prior to initial use as well as each additional use.
- Only use endoscopes which are in perfect condition.
- Clean / disinfect and sterilise the endoscope prior to initial use as well as each additional use of the endoscope. Contaminants on the irradiation surface of the illumination fibres can burn in during use, which impacts image quality.
- Ensure that the proximal end of the endoscope is dry to prevent the endoscope from fogging up during the examination / procedure.
- Ensure that no parts are missing or loose.
- Ensure that there are no residual cleaning agents or disinfectants on the endoscope.
- Inspect the entire endoscope, particularly the sheath, for contaminants and damage of any type, such as dents, scratches, cracks, bending and sharp edges.
- Inspect distal end, proximal end and irradiation surface of the illumination fibres for contamination and scratches. Make contaminants and scratches visible using light reflexes. Hold the connection of the optical fibres against the light and inspect whether the optical fibres illuminate evenly at the distal end.
- Check image quality: The image may not be blurry, clouded or dark. If necessary, remove deposits on the optical end surface using polishing paste provided.
- For endoscopes with locking device: Inspect between the sheath and the main part for contaminants and damage to ensure a fixed and secure connection.
- For endoscopes with working channel: Inspect all components of the adapter with stopcock for function and damage.
- Inspect free passage of the working channels.

7.2 Provisioning

- Clean / disinfect and sterilise the endoscope prior to initial use as well as each additional use of the endoscope, see "Processing" on page 19.
- Ensure that the proximal end of the endoscope is dry to prevent the endoscope from fogging up during the examination / procedure.
- If necessary mount adapter for illumination fibre, see "Assembly" on page 23.
- Mount illumination fibre (see manufacturer's specifications).
- If required, adapt the camera (see manufacturer's specifications).

8. Use

WARNING

Risk of burns

The optical fibres emit high-energy light at the distal end of the endoscope. This can cause the temperature of the body tissue to rise to 41 °C.

- Avoid direct contact of the distal end with body tissue or flammable materials as it can cause burns.
- Reduce the light intensity of the cold light source when working near body tissue or flammable materials.

For endoscopes with working channel:

- If an instrument is inserted, do not move the lever on the stopcock because otherwise the instrument could be damaged or sheared off.
- Prepare the endoscope for processing immediately after use to prevent surface drying of the contaminants.

9. Processing

9.1 Warnings and precautions

- Observe the country-specific regulations and laws concerning cleaning of endoscopes for patients with Creutzfeld-Jakob-Disease, CJK-on-spec or its possible variants, Bovine Spongiform Encephalopathy or Transmissible Spongiform Encephalopathy.

Limitations on reprocessing

- Using multiple sterilization methods may significantly reduce the performance of the endoscope and is not recommended.
- Damage incurred by improper processing will not be covered by the warranty.

9.2 Cleaning and disinfection

Transport

- Always store endoscope securely and transport it to processing in a closed container to prevent damage to the endoscope and the contamination of the environment.

Manual cleaning / pre-cleaning and chemical disinfection

- The manual cleaning procedure must be performed prior to the automated cleaning.

- Do not use fixating cleaning agents or hot water (> 40 °C) as it can cause fixation of the contaminants and jeopardise successful cleaning.
- Do not scratch contaminants off with hard objects as this can cause damage to the endoscope.

- For endoscopes with stopcocks: All components of the stopcocks are dismantled, see chapter 10.
- Remove coarse contamination from the endoscope. Clean the endoscope with a soft brush under cold tap water until all visible contaminants have been removed.
- For working channels, holes and threads. Rinse with a water pistol for at least 10 seconds at a pressure of at least 3.8 bar (absolute).
- Disinfect the endoscope. In doing so, observe the specifications of the disinfectant solution manufacturer regarding temperature, concentration and application time.

Non-compliance with the manufacturer's specifications can result in damage to the endoscope.

- Rinse endoscope with deionised water.
- Dry working channel with compressed air.
- Dry endoscope with lint free cloths.
- Carry out visual inspection, function check and servicing.

EN

Machine cleaning and thermal disinfection

The endoscopes are suitable for prevalent machine methods of cleaning and thermal disinfection. In doing so, use gentle cycles for endoscopes and suitable cleaning agents and disinfectants. The instructions of the machine, cleaning agent and disinfectant manufacturers must be observed. The cleaning and disinfectant result must be confirmed by the machine, cleaning agent and disinfectant manufacturers in cooperation with the user.

- Fix the endoscopes to the corresponding tray in such a way that a successful cleaning procedure is possible and that damage is prevented during cleaning, see pictures of the trays in chapter 9.6.
- Place the tray in an automated washer.
(The automated cleaning was validated with the cleaning unit Miele Professional G 7836 CD. For use with another washer/disinfector make sure that it is also in compliance with EN ISO 15883-2.)
- Program the washer with the following parameters and activate the wash.

Step	Process step	Reagents	Time	Temperature
1	Pre-cleaning	Tap water	3 min	Cold tap water ($\leq 23\text{ }^{\circ}\text{C}$)
2	Drain	—	—	—
3	Cleaning	Tap water • Dosing: 0.5% neodisher MediClean forte at $55\text{ }^{\circ}\text{C}$	5 min	$55\text{ }^{\circ}\text{C}$ (set point)
4	Drain	—	—	—
5	Rinsing	Deionised water	1 min	Cold deionised water ($\leq 23\text{ }^{\circ}\text{C}$)
6	Drain	—	—	—
7	Rinsing	Deionised water	1 min	Cold deionised water ($\leq 23\text{ }^{\circ}\text{C}$)

Disinfection

Thermal disinfection has been validated using the following parameters:

- Time: 95 seconds
- Temperature: $95\text{ }^{\circ}\text{C}$
- Ensure that the exteriors of the endoscope are dry. If necessary, dry with lint free cloths.
- Carry out visual inspection, function check and servicing.

9.3 Sterilisation

i Prior to each sterilisation, the endoscopes must be cleaned and disinfected according to the methods in these Instructions for use.

- Double wrap the tray prior to sterilisation according to AAMI ST 79. Sterilise endoscopes in suitable packaging to prevent subsequent contamination.
(The validation was performed with Broemeda Amcro Flexibles sterilisation wraps $900 \times 900\text{ mm}$ (Broemeda item number 68170990).)

Steam sterilisation (autoclaving)

In general, users are responsible for the validation of their processes.

- When selecting the processing method, observe the valid national hygienic regulations and local provisions for hospital hygiene.

Comply with specified process parameters. The specified parameters have been validated to ensure the sterility of the endoscopes. Deviating process parameters could damage the endoscope. In this case, the guarantee and warranty shall become void.

- Sterilise endoscopes.
- When the sterilisation process has ended, allow the endoscopes to cool gradually to room temperature.

Fractionated pre-vacuum method ¹⁾

The following process has been validated in accordance with ISO 17665:

Cycle type	Temperature	Exposure time ²⁾	Minimum dry time	Minimum cool time
UK prevacuum/ pulsating vacuum	134 °C / 273 °F	3 min	30 min	30 min
Prevacuum/ pulsating vacuum	132 °C / 270 °F	4 min		
Prevacuum/ pulsating vacuum	134 °C / 273 °F	18 min		

- The sterilisation validation was validated with the steam steriliser Selectomat HP 666-1HR of MMM.
- 18 minutes is the maximum exposure time that should be used for the prevacuum cycle.

Packaging

The products are delivered non-sterile in sealed plastic or in a protective box/foam packaging. Transport packaging is not suitable for sterilisation.

- Pack the products into suitable sterilisation packaging systems (e.g. STERICLIN pouch used during sterilisation validation) acc. to ISO 11807 and/or AAMI / ANSI ST77.2006 in order to be sterilised.

9.4 Special precautions: Pathogens of Transmissible Spongiform Encephalopathy

A comprehensive explanation of the necessary preventative measures with regard to agents of Transmissible Spongiform Encephalopathy (TSE) would go beyond the scope of this document.

It is assumed that pathogens of the Creutzfeldt Jakob Disease cannot be killed using normal disinfection and sterilisation processes. Therefore, the standard methods for decontamination and sterilisation are not sufficient if there is a risk of transferring Creutzfeldt Jakob Disease.

In general, only tissue with a low potential of TSE infection comes into contact with surgical instruments. In spite of this, special preventative measures must be taken for instruments which are used to treat patients with a known or suspected infection of TSE, as well as for patients at risk.

9.5 Processing restrictions

Repeated processing has only minimal effect on the endoscopes. The lifetime of the units is usually determined by wear and damage.

The endoscope can be damaged if the manufacturer's specifications are not observed.


- Do not clean endoscope in an ultrasonic bath.

EN


9.6 Tray configuration

**Tray configuration 4 mm set automated cleaning
(2.7 mm set similar)**


- 1 Adapters, stopcock plug and nut
- 2 8300036854
- 3 8300036862
- 4 84-20-5002
- 5 8300036857
- 6 Stopcock bridge

**Tray configuration 4 mm set steam sterilisation
(2.7 mm set similar)**


- 1 8300036854
- 2 8300036862
- 3 84-20-5002
- 4 8300036857

Tray configuration 1.9 mm set automated cleaning


- 1 Adapters, stopcock plug and nut
- 2 8300036850
- 3 8300036860
- 4 84-20-5002
- 5 8300036859
- 6 Stopcock bridge

Tray configuration 1.9 mm set steam sterilisation

- 1 8300036850
- 2 8300036860
- 3 84-20-5002
- 4 8300036859

10. Assembly


- 1 Connection for illumination fibre, type ACMI
- 2 Adapter type Wolf
- 3 Adapter type Storz / Olympus

- If necessary mount adapter for illumination fibre
- Ensure that the irradiation surface of the illumination fibre is clean.
- Mount illumination fibre (see manufacturer's specifications).
- If required, adapt the camera (see manufacturer's specifications).

EN

For endoscopes with working channel:

- a Union nut
- b Stopcock
- c Sealing cap, green
- d Stopcock nut

i In order to ensure sterility, only use grease which is suitable for medical instruments for the stopcock.

- Lubricate stopcock b.
- Mount stopcock b and fix with stopcock nut d.
- Remove excess grease.

11. Disassembly

WARNING

Risk of burns!

Allow the illumination fibre to cool sufficiently before removing it. The ends can become very hot and cause serious burns.

- Remove illumination fibre.

- Do not remove the ocular funnel because otherwise the endoscope will be damaged.

- Unscrew existing adapters.
- For endoscopes with working channel, disassemble the adapter with stopcock.
 - Remove sealing cap c.
 - Loosen union nut a.
 - Loosen stopcock nut d.
 - Disassemble stopcock b.

12. Storage

Unsterile metal units must be stored in a clean, dry environment. The storage time of unsterile units is not limited; the units are made of a non-degradable material which maintains its stability when stored under the recommended conditions.

As long as endoscopes are stored unsterile in the original packaging, the following storage conditions apply:

- Avoid direct sunlight.
- Store endoscope either in the original packaging or in a screen tray/container.
- Ensure that the endoscope is stored securely.
- Observe the respective valid national provisions when storing in a sterile condition.

13. Ambient conditions

	Operation	Transport and storage
Temperature	10 °C - 40 °C	-40 °C - 70 °C
Relative humidity	5 % - 75 %	5 % - 95 %
Atmospheric pressure	70 kPa - 106 kPa	70 kPa - 106 kPa

14. Service and maintenance

Henke-Sass, Wolf does not supply original parts to independent workshops or other endoscope manufacturers. Thus only Henke-Sass, Wolf is in a position to carry out repairs using original parts. The original technical specifications and the operational safety of the endoscope can only be guaranteed by using original parts. The warranty for Henke-Sass, Wolf products shall become void if repairs are carried out by an workshop not authorised by Henke-Sass, Wolf. In this case Henke-Sass, Wolf is also no longer responsible for the technical specifications or safety of the product.

- Have the endoscopes repaired by Henke-Sass, Wolf only. To this purpose send the defective endoscope to the address of the Service Centre of Zimmer Biomet.
- Clean, disinfect and sterilise the endoscope thoroughly prior to returning it for repair.
- Ideally, send in the endoscope in its original packaging. If this is not possible, package the endoscope to secure it for transport.

Henke-Sass, Wolf is not liable for damage resulting from improper shipping.

15. Disposal

- Observe country-specific regulations and laws for the disposal of medical products.

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© copyright Henke-Sass, Wolf