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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of REN FENG of MD PACIFIC (TIANJIN) BIOTECHNOLOGY CO., LTD. and the seal of the said company on the annexed DOCUMENT are genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Wang Shuping

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(Date: Jun. 28, 2020)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

2
APPROVED BY

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
User manual

Semi-automatic coagulometer "Technology Solution 4»

Please read the operating instructions carefully before installation and use.

Version: 1.0. Limited Liability Company "Laboratory Solutions". 656037, Russia, Barnaul,
Kalinina ave., 116/10, office 108, т. 8-800-551-61-71, sales@lab-solution.ru

How to use this manual

The semi-automatic coagulometer "Technology Solution 4" must be operated by trained operators.

Please read the "User Manual" carefully before using it!

Thank you for being the owner of the semi-automatic coagulometer "Technology Solution 4"! Please read the manual carefully before using the device. To achieve optimal results, it is necessary to become familiar with the principle of operation of the equipment before starting clinical trials. This guide will provide information about installation, daily use, maintenance, and other aspects of working with the semi-automatic coagulometer "Technology Solution 4". Keep all packaging materials – they are needed for storage, transportation, and maintenance of the product. For all technical questions, please contact the distributor.

Warning: if you ignore the safety precautions when working with the device, this can lead to distortion of the analysis results, damage to the analyzer parts, and even a threat to the operator's health.

Note: The device is used for clinical research. A clinical diagnosis based on the results of the analysis should be made by a doctor in accordance with the clinical symptoms and a number of other results.

Notification

Limited Liability Company "Laboratory Solutions" reserves the last right to amend the operating manual.

Illustrations are not intended for other purposes. Without written permission of LLC "Laboratory Solutions" it is prohibited to reproduce or modify this manual.

The manufacturer is responsible for the safety, reliability and stability of the product only subject to the following requirements:

- assembly, re-adjustment, rework and repair shall be carried out by persons with special permission of LLC "Laboratory Solutions";
- the product must be operated in accordance with this manual;
- all electrical devices connected to the analyzer must meet the requirements of international standards.

Note

Specialists in the field of laboratory diagnostics are allowed to work with the device:

laboratory managers, QLD doctors , laboratory technicians, medical laboratory technicians (medical assistant), laboratory technicians, medical technologists.

Attention

If the user is unable to provide proper operation/maintenance, the device may fail and even pose a threat to the user's health.

Ensure that the device is operate under the conditions specified in this manual. Otherwise, the device may malfunction, distortion of the analysis results, damage to the analyzer parts, and injury.

Symbols on the package:

Manufacturer

Fragile freight

There is a beating load in the transport package, so it should be handled with special care

Top

The transport package shall be on the specified side upwards

Protect from rain

The packaged product is sensitive to moisture

Serial number

Medical device for diagnostics *in vitro*

Symbols on marking:

	Serial number
---	---------------


	Date of manufacture
	Manufacturer
	Medical device for diagnostics <i>in vitro</i>
	Biological hazard

Symbols in the analyzer panel:

O (OFF)	To turn off the power
I (ON)	To turn on the power
a. s. 220 V / 50 Hz	The parameters of the power supply
FUSE	Fuse icon
0.5A×2	Nominal value and number of fuses The analyzer has 2 fuses with a diameter of 5 mm, a length of 20 mm and a current of 0.5 A.

Marking layouts for the semi-automatic coagulometer "Technology Solution 4»

Layouts of packaging marking for semi-automatic coagulometer "Technology Solution 4»

Safety precautions and recommendations

The semi-automatic coagulometer "Technology Solution 4" is design for in vitro diagnostics. Human blood plasma stabilized with an anticoagulant (sodium citrate) is use for analysis. Any other material is considered unintended. These operating conditions also provide for compliance with the maintenance procedures described in the instructions.

Content

How to use this manual	2
Symbols on the package	3
Symbols on marking	3
Symbols in the analyzer panel	4
Marking layouts for the semi-automatic coagulometer "Technology Solution 4»	4
Layouts of packaging marking for semi-automatic coagulometer "Technology Solution 4»	5
Safety precautions and recommendations	5
1. Introduction	9
1.1 Features of the device	9
1.2 The structure and components of the analyzer	9
1.3 Scope of the analyzer	10
2. Parameters and features	11
2.1 Technical specifications	11
2.1.1 Electromagnetic radiation	12
2.1.2 Electromagnetic noise stability	13
2.2. Description of the composition of the semi-automatic coagulometer "Technology Solution 4»	16
3. Configurations and functions	23
3.1 Configuration	23
3.1.1 Front panel	23
3.1.2 Rear panel	24
3.2 Specifications	24
4. Functions	25
4.1 Functional scheme of the device	25
4.2 Function buttons	26
5. Operation of the device	27
5.1 Installation	27
5.2 Turning on	27
5.3 Beginning of use	27

5.4 Get zero point values	27
5.5 Choice of the test	28
5.5.1 Choice of one test	28
5.5.2 Choice of several tests	29
5.5.2.1 Editing group of tests	30
5.5.2.2 Choice of group of tests	31
5.6 Preparation for measurement	32
5.6.1 Changing the ID number	32
5.7 Preliminary incubation of samples	33
5.7.1 Beginning of a preliminary incubation	33
5.8 Test	34
5.8.1 The beginning of the measurement	34
5.8.2 End of the test	35
5.8.2.1 Normal end of the test	35
5.8.2.2 Abnormal end of the test	37
5.9 Main menu	37
5.10 Management of reports	38
5.10.1 Printing of reports	39
5.10.2 Deleting a single report	41
5.10.3 Deleting all reports	42
5.10.4 Uploading reports	44
5.11 Parameters	45
5.11.1 Choice of parameter	45
5.11.2 Setting of parameters	46
5.12 System settings	48
5.12.1 Date and time	48
5.12.2 Choice of language	49
5.12.3 Print settings	50
5.13 Calibration	51
5.13.1 Set the calibration curve for the FIB test	52
5.13.2 Set the calibration curve for the PT test	54

5.14 Quality control	57
5.14.1 Representation of quality control on a graph	58
5.15 Service menu	59
5.16 Technical support	59
5.17 Printer	59
5.17.1 Installing paper in the printer	60
5.17.2 Removing paper from the printer	60
5.17.3 Thermal printer maintenance	60
5.18 Connection with PC	60
6. Work Precautions	61
7. Important notes and maintenance	62
7.1 Important notes	62
7.2 Maintenance	63
8. System errors and their elimination	64
9. Manufacturer's warranty	65
10. Packaging, storage and transportation	66

1. Introduction

1.1 Features of the device

Medical device semi-automatic coagulometer "Technology Solution 4" is a semi-automatic coagulometer with a mechanical detection method (double magnetic chain with a ball) and is used for analyzing the time and activity of blood clotting factors, analyzing the activity of anticoagulant and fibrinolytic systems in the field of clinical diagnostics of diseases associated with hemostasis pathology, as well as for monitoring thrombolytic and anticoagulant therapy. The analyzer is easy to use, stable performance, and accurate measurements. The principle of operation of this coagulation analyzer is a double magnetic circuit and a magnetic ball that fixes the formation of a clot when its vibrations fade. The test result is not affected by hemolysis, chylousness, icteritiousness, the turbidity of the plasma and other factors.

The analyzer uses a clear interface, each step of the operation is described in this manual. Intelligent temperature control and time control provide a simple and accurate way to get results. Built-in thermal printer simplifies user experience. An optimized user operating system provides a simple testing process. Pre-installed data transfer systems allow you to connect to the Laboratory Information System (LIS) through the RS232 port.

The product includes:

1. Semi-automatic coagulometer "Technology Solution 4" - 1 pc.
2. User manual - 1 pc.
3. The passport of the device - 1 pc.
4. Power cord - 1 pc.
5. Fuse - 2 pcs.
6. Dispenser for magnetic balls - 1 pc.
7. Paper for thermal printer - 2 rolls.
8. Stylus - 1 pc.
9. Automatic pipette (if necessary) - 1 pc.

1.2 The structure and components of the analyzer: display, thermal printer, operating system and measurement system.

To obtain correct analysis results and correct operation of the analyzer, only single-use specialized cells with the following characteristics must be used:

- Volume of 500 μ l;

- Total height of the cuvette- 29.85 ± 0.15 mm;
- Overall size (length×width) of the landing part of the cuvette - $10.80 \pm 0.1 \times 7.00 \pm 0.1$ mm;
- Internal size (length×width) of the bottom of the cuvette - $9.65 \pm 0.1 \times 4.50 \pm 0.1$ mm;
- External overall size (length×width) of the bottom of the cuvette - $10.55 \pm 0.1 \times 4.50 \pm 0.1$ mm;
- The distance of the measuring area from the stop plane is 3.25 ± 0.1 mm;
- External overall size (length×width) of the measuring part of the cuvette - $10.60 \pm 0.1 \times 4.60 \pm 0.1$ mm;
- The wall thickness in the measuring part of the cell is not more than 0.45 ± 0.05 mm;
- The radius of the bottom forming part in the measuring part of the cuvette is 11.32 ± 0.1 mm;
- The distance from the top plane to the bottom of the cell at the extreme point is 26.60 ± 0.1 mm;
- The distance from the upper plane to the bottom of the cell at the midpoint is 28.80 ± 0.1 mm •
- Single-use measuring ball.

1.3 Scope of the analyzer

The test parameters measured by the device are shown in the table:

Table 1

Name	Measured parameter	Method of analysis	Units of measure
Prothrombin time	PT	Coagulation	sec (second)
Activated partial thromboplastin time	APTT	Coagulation	sec (second)
Fibrinogen (modified Clauss method)	FIB	Coagulation	g / L (g / l)
Thrombin time	TT	Coagulation	sec (second)
Lupus anticoagulant LA +	LA+	Coagulation	sec (second)
Lupus anticoagulant LA-	LA-	Coagulation	sec (second)
Protein C PARUS1	PARUS1	Coagulation	sec (second)
Protein C PARUS2	PARUS2	Coagulation	sec (second)
Coagulation Factor VIII	F8	Coagulation	%
Coagulation Factor IX	F9	Coagulation	%

Coagulation Factor X	F10	Coagulation	%
Coagulation Factor XI	F11	Coagulation	%
Coagulation Factor XII	F12	Coagulation	%

2. Parameters and features

2.1 Technical specifications

- Normal operating temperature range: 10 °C ~ 30 °C
- Atmospheric pressure: 86 kPa ~ 106 kPa
- Voltage: 220V
- Frequency: 50/60 Hz
- Relative humidity: ≤70%
- Input power: 70 W
- Preheating time: no more than 30 minutes (at room temperature from 20 °C to 37 °C)
- Temperature control: incubation and measuring cells maintain a temperature of 37.0 °C ± 1.0 °C.
- Reagent cell heating is 37.0 °C ± 1.0 °C.
- The discrepancy in the channels: the maximum discrepancy in the test results on different channels does not exceed ≤5%
- Measurement accuracy: the relative bias of the FIB results is less than ± 10.0%.
- Reproducibility: refer to the following information.

Table 2

Test	The coefficient of variation CV	
	Normal sample	Pathological sample
PT (sec)	≤ 5.0%	≤ 10.0%
APTT (sec)	≤ 5.0%	≤ 10.0%

FIB (g/l)	≤ 10.0%	≤ 10.0%
TT (sec)	≤ 10.0%	≤ 10.0%
VIII (%)	≤ 10.0%	≤ 10.0%
IX (%)	≤ 10.0%	≤ 10.0%
X (%)	≤ 10.0%	≤ 10.0%
XI (%)	≤ 10.0%	≤ 10.0%
XII (%)	≤ 10.0%	≤ 10.0%

- Continuous working time: at least 24 hours.
- Dimensions : 420×345×155 mm
- Security class for software according to EN 62304 – A
- Electrical safety according to EN 61010-1:
- Overvoltage category I.
- Degree of contamination 1.
- Protection against electric shock in a medical device is provided by protective earthing.
- Group of dependence on perceived mechanical effects according to GOST R 50444-92: 2-wearable, portable and mobile, not intended for use when carrying and moving within a medical institution
- The level of protection of the analyzer against dangerous penetration of water or solid particles-IP00
- No laser radiation – mechanical detection method.

2.1.1 Electromagnetic radiation

Table 3

The semi-automatic coagulometer "Technology Solution 4" is intend for use in the electromagnetic environment defined below. The user must ensure that the device is operate under these conditions.		
Electromagnetic emission test	Accordance	Electromagnetic environment - instructions
Radio interference according to CISPR 11	Group 1	The Technology Solution 4 semi-automatic coagulometer uses radio frequency energy only to perform internal functions.

		The level of radio frequency interference emissions is low and probably will not lead to malfunctions of nearby electronic equipment
Radio interference according to CISPR 11	Class A	Semi-automatic coagulometer "Technology Solution 4" belongs to class A (class A equipment): Equipment intended for use in locations that are not related to residential areas, as well as in locations where the equipment is not directly connected to low-voltage distribution networks that supply energy to residential buildings.
Harmonic components of current consumption in accordance with IEC 61000-3-2	Do not use	
Voltage fluctuations and flicker according to IEC 61000-3-3	Do not use	

2.1.2 Electromagnetic noise stability

Table 4

Test for noise stability	IEC 60601 test level	Electromagnetic environment - instructions
Electrostatic discharges (ESD) according to IEC 61000-4-2	$\pm 2 \text{ kV}$, $\pm 4 \text{ kV}$, $\pm 6 \text{ kV}$ - contact discharge $\pm 2 \text{ kV}$, $\pm 4 \text{ kV}$, $\pm 8 \text{ kV}$ - air discharge	The floors of the room should be wooden, concrete or tiled. If the floors are covered with synthetic material, then the relative humidity should be at least 30%
Nanosecond impulse noise according to IEC 61000-4-4	$\pm 1 \text{ kV}$ - for power lines $\pm 0.5 \text{ kV}$ - for input-output lines	The quality of electrical energy in the building's electrical network must meet the typical conditions of a commercial or hospital setting

High energy microsecond impulse noise according to IEC 61000-4-5	± 1 kV - when applying wire-to-wire interference ± 2 kV - when applying wire-to-ground interference	The quality of electrical energy in the building's electrical network must meet the typical conditions of a commercial or hospital setting
Voltage dips, interruptions and changes in input power lines according to IEC 61000-4-11	0% U_N (voltage drop 100%) within 1 period 40% U_N (voltage drop 60% U_N) for 5 periods 70% U_N (voltage drop 30% U_N) for 25 periods <5% U_N (voltage interruptions >95% U_N) for 250 periods	The quality of electrical energy in the building's electrical network must meet the typical conditions of a commercial or hospital setting.
Industrial frequency magnetic field (50 Hz) according to IEC 61000-4-8	3 A/m	Industrial frequency magnetic fields must be measured at the designated installation site to ensure that the field strength is low enough
Conducted interference induced by radio frequency electromagnetic fields in accordance with IEC 61000-4-6	3 V (RMS value)	The recommended spatial diversity is: $d = 1,2\sqrt{P}$
Radiated radio frequency electromagnetic field according to IEC 61000-4-3	3 V / m in the band from 80 MHz to 2.5 GHz 1 V / m in the band from 2.0 GHz to 2.7 GHz	$d = 1,2\sqrt{P}$ (from 80 to 1000 MHz); $d = 1,2\sqrt{P}$ (from 80 to 1000 MHz); $d = 2,3\sqrt{P}$ (from 2.0 to 2.7 GHz).

--	--	--

- Notes:** 1. At frequencies of 80 and 1000 MHz, a higher value of the field strength is used.
- 2 The above expressions are not applicable in all cases. Electromagnetic propagation is affected by absorption or reflection from structures, objects, and people.

Table 5

Recommended values of spatial separation between portable and mobile radio frequency communication devices and semi-Automatic coagulometer " Technology Solution 4»			
The semi-automatic coagulometer "Technology Solution 4" is intend for use in an electromagnetic environment where the levels of radiated interference are monitored. The user of the semi-Automatic coagulometer "Technology Solution 4" can avoid the influence of electromagnetic interference by ensuring a minimum spatial separation between portable and mobile radio frequency communication means (transmitters) and the device, as recommended below, taking into account the maximum output power of the communication means.			
Nominal maximum output power of the transmitter, W	Spatial variation, m, depending on transmitter frequency		
	$d = 1,2\sqrt{P}$ in the band from 80 MHz to 1000 MHz	$d = 1,2\sqrt{P}$ in the band from 1.4 to 2.0 GHz	$d = 2,3\sqrt{P}$ in the band from 2.0 GHz to 2.7 GHz
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
Notes: 1 At frequencies of 80 and 1000 MHz, a higher value of the field strength is applied. 2 These expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people. 3 When determining the recommended spatial spacing values for transmitters with a nominal maximum output power not specified in the table, the following expressions are substitute for the nominal maximum output power in watts specified in the manufacturer's documentation.			

2.2. Description of the composition of the semi-automatic coagulometer "Technology Solution 4»

Table 6

No. p/ p	Name	Appointment	Technical parameters
1.	Power cord	Provides electricity to the device	Length 1740 mm $\pm$ 10 mm. Voltage: up to 500 V. Diameter: 8 mm $\pm$ 1 mm.
2.	Fuse	It performs a protective function, protects the electric circuit and its elements from overheating and fire when a high current flows.	Diameter: 5mm Length: 20 mm. Current strength: 0.5 A. Voltage: 220 V.
3.	Dispenser for magnetic balls	Designed for dispensing balls into measuring ditches	Length: 90 mm $\pm$ 1 mm Diameter: 47 mm $\pm$ 1 mm
4.	Thermal Printer Paper	Designed to display test results on paper	Maximum tape diameter: 40 mm $\pm$ 1 mm Tape width: 57 mm. Paper weight 80 g / m2.
5.	Stylus	Used to control the touch interface of the coagulometer	Type: Capacitive Length: 100 mm $\pm$ 1 mm.
6.	Automatic pipette	It used to add a plasma sample and reagent, and start testing.	Length: 222 mm $\pm$ 1 mm Dosing Range: 20-200 μ l

Functional characteristics of the semi-Automatic coagulometer "Technology Solution 4" in accordance with GOST R ISO 18113-3-2015:

- Throughput of about 80 tests per hour (prothrombin time (RT));
- Transfer is excluded. Separate disposable tips for dispensers and disposable cuvettes are used;
- Cross-contamination is excluded. Separate disposable tips for dispensers and

disposable cuvettes are used;

- Sample volume from 20 to 200 μl , depending on the method;
- Reagent volume from 20 to 200 μl , depending on the method;
- Measurement time from 5 to 200 seconds;
- Measurement temperature 37 degrees Celsius;
- There is no linear absorption of the photometer – mechanical detection method;
- No wavelength-mechanical detection method.

Information about the classification of Automatic pipette in accordance with GOST 28311-89:

Depending on:

- degree of mechanization of the process of forming the dose volume-with manual drive;
- the number of simultaneously generated quantities of doses - single channel;
- possibility of changing the dose volume - with a variable volume;
- locations of the dispenser during operation-held in the hand of the laboratory assistant;
- functions performed - for reagents;
- contact with the liquid metering unit - with an air gap;

The main parameters of the Automatic pipette in accordance with GOST 28311-89:

- The nominal value of the highest dosage limit for a dispenser with variable dose volumes is $2 * 10^2 \mu\text{l}$.
- The nominal value of the lowest dosage limit for a dispenser with variable dose volumes of no more than 10 %.
- Discreteness of setting the dose volume for dispensers with a variable dose volume of $5 * 10^0 \mu\text{l}$.
- The limit of the characteristics of deviations in the volume of doses set for the dispenser when dispensing distilled water does not exceed the following values:
- The limit of the permissible relative deviation of the arithmetic mean of the actual volume of the dose from the nominal, % (dispenser with an air gap) - $\pm 1.5\%$.
- The limit of the permissible relative average square deviation of the actual dose volume with a confidence probability, % (air-gap dispenser) - 2%.

Technical specifications of Automatic pipettes in accordance with GOST 28311-89:

Characteristics of deviations of volumes of doses normalized by setting:

the limits (positive and negative) deviation of the arithmetic mean of the actual dose volume from the nominal and maximum permissible standard deviation of the actual dose volume with a confidence level selected from a number of 0,90; 0,95; 0,98; 0,99.

Nominal inspected volume, μl	Limit of permissible relative deviation of the arithmetic mean value of the actual dose volume from the nominal one, $\sigma\%$	Limit of allowable relative mean quadratic deviation of actual dose volume, $\Delta\%$
20	$\pm 2,0$	$\pm 3,0$
25	$\pm 2,0$	$\pm 3,0$
30	$\pm 2,0$	$\pm 3,0$
50	$\pm 2,0$	$\pm 2,5$
100	$\pm 1,5$	$\pm 2,0$

Number of channels	Range, (μl)	Volume, (μl)	Accuracy, (%)	Reproducibility, (%)
1	20-200	200	0,6	0,15
		100	0,8	0,30
		20	2,0	0,80

The change in the arithmetic mean value of the actual dose volume when the influencing value changes under operating conditions or the liquid being dosed does not exceed 0.35 of the set limit.

- The limits of variation of deviations of the actual dose value when resetting the dose volume do not exceed 0.5 absolute value of the deviation limit of the arithmetic mean of the actual deviation of the dose volume from the nominal one.

Reliability requirements

- Reliability requirements must be set in the manufacturer's technical documentation for specific types of dispensers in the form of standardized values of reliability, durability and maintainability indicators in accordance with GOST 23256-86.
- Trouble-free – no less than the operating time during the warranty period.

- Durability (product service life) - 5 years.
- The product is not repairable.
- Depending on the possible consequences of failure during operation, the dispenser belongs to class B according to GOST 23256-86.
- The dispenser and its parts, held in the hand of the laboratory assistant, have a mass of no more than 0.35 kg.

Requirements for resistance to external influences

- The dispenser is resistant to mechanical factors for product group 2 according to GOST 20790.
- The dispenser during operation is resistant to climatic factors according to GOST 20790 for the type of climatic version of UHL 4.2 according to GOST 15150.
- The metering device is resistant to climatic factors during transportation: the metering device of climatic version UHL 4.2-for group of transportation conditions 1;
- The dispenser is resistant to climatic factors during storage: the dispenser of climatic version UHL 4.2 - for group of storage conditions 2;
- The outer surfaces of the dispenser are resistant to disinfection and, if necessary, to sterilization.
- The materials of the dispenser parts (or their coatings) that come into contact with the dispensed liquids are resistant to their effects.
- For a single-channel dispenser with a dose volume of up to 1 ml held in the hand of a laboratory assistant, the greatest effort required to form the dose volume using the fingers does not exceed 25 N.
- Aesthetic indicators of the dispenser ensure informational expressiveness, integrity of the composition, and perfection of production performance.

Functional characteristics of the Automatic pipette in accordance with GOST R ISO 18113-3-2015

- Throughput is determined by the performance of the device;
- Transfer is excluded. Separate disposable tips for dispensers are used;
- Cross-contamination is excluded. Separate disposable tips for dispensers are used;
- Sample volume from 20 to 200 μl , depending on the method;
- Reagent volume from 20 to 200 μl , depending on the method;
- The time of measurement, the pipette is not held;
- Temperature measurement – pipette is not carried out;

- Linearity of absorption – pipette is not carried out;
- The wavelength is not pipetted.

Kits and reagents are used to operate the medical device as intended

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№	Article	Name	Registration certificate
1.	652	Reagent kit for determining activated partial thromboplastin time (APTV-El-TEST) according to TU 9398-290-42349142-2009	ФСР 2010/06655 от 22.06.2010
2.	731	Reagent kit for determining activated partial thromboplastin time (APTV-El-TEST) according to TU 9398-057-42349142-2015	РЗН 2016/4444 от 08.07.2016
3.	748, 735	A set of reagents for determining prothrombin time (Techplastin test) according to TU 9398-056-42349142-2015	РЗН 2016/4433 от 08.07.2016
4.	131, 607, 140, 608	A set of reagents for determining prothrombin time (Techplastin test) according to TU 9398-289-42349142-2005	ФСР 2010/08020 от 22.06.2010
5.	711	A set of reagents for determining the concentration of fibrinogen on automatic and semi-automatic coagulometers (MultiTech-Fibrinogen) according to TU 9398-029-42349142-2012	ФСР 2012/13639 от 29.06.2012
6.	742	Reagent for the study of hemostasis (Thrombin) according to TU 9398-060-42349142-2015	РЗН 2016/4521 от 01.08.2016

7.	151, 609, 610	Set of reagents for determination of thrombin time (Thrombo test) according to TU 9398-230-42349142-2003	ΦCP 2007/015729 от 22.06.2010
8.	193	A set of reagents for determining lupus anticoagulant (Express-Lupus test) according to TU 9398-027-42349142-2012	ΦCP 2012/13467 от 24.05.2012
9.	164	A set of reagents for determining disorders in the protein C system (PARUS-TEST) according to TU 9398-291-42349142-2009	ΦCP 2010/06656 от 22.06.2010
10.	274	A set of reagents for determining the activity of factor VIII in blood plasma (Tech-Factor VIII-test) according to TU 9398-281-42349142-2004	ΦCP 2008/02322 от 22.06.2010
11.	679	A set of reagents for determining the activity of factor IX in blood plasma (Tech-Factor IX-test) according to TU 9398-017-42349142-2011	ΦCP 2011/11955 от 22.09.2011
12.	027	Tris-HCl Buffer according to TU 9398-046-42349142-2015	P3H 2016/4445 от 08.07.2016
13.	OTXY13	In vitro reagents for hemostasis research. 34. Factor X-deficient plasma.	ΦC3 2009/03943 от 04.02.2020
14.	OSDF13	In vitro reagents for hemostasis research. 33. Factor IX-deficient plasma.	ΦC3 2009/03943 от 04.02.2020


15.	OSDG13	In vitro reagents for hemostasis research. 29. Factor XII-deficient plasma.	ΦC3 2009/03943 от 04.02.2020
16.	773	Calibration Plasma (Multitech Calibrator) according to TU 9398-064-42349142-2016	P3H 2017/5432 от 27.02.2017
17.	774	Control plasma with a normal range of values (Technoclot N) according to TU 9398-065-42349142-2016	P3H 2017/5653 от 18.04.2017
18.	775	Control plasma with a pathological range of values (Tekhnoklot P) according to TU 9398-066-42349142-2016	P3H 2017/5640 от 12.04.2017
19.	776	Control plasma with a normal range of values (Tech control N) according to TU 9398-067-42349142-2016	P3H 2017/5723 от 05.05.2017
20.	777	Control plasma with a pathological range of values (Tech control P) according to TU 9398-068-42349142-2016	P3H 2017/5720 от 05.05.2017

3. Configurations and functions

3.1 Configuration

3.1.1 Front panel

Structure of the front panel of the semi-automatic coagulometer "Technology Solution 4»

Pic. 3-1 Configuration of the front panel of the analyzer "Technology Solution 4»


- Built-in thermal printer and switch: press the button to open the cover.
- LCD touch screen: 7 inch TFT touch screen.
- Automatic pipette: used to add reagent and start testing.
- Reagent position: 2 reagent positions without heating.
- Heated reagent positions: 2 reagent positions heated to a temperature of 37.0 °C.
- Measuring channels: 4 measuring channels from top to bottom from 1st to 4th.
- Incubation cells: 16 (4x4) incubation cells with a constant temperature of 37.0 °C.

The device case can be produce in two color versions: white and gray.

3.1.2 Rear panel

Pic. 3-2 shows the rear panel of the analyzer.

- Switch: when the switch is in position "I", the device is turned on; when the switch is in the "O" position, the device is turned off.
- Power connector: connect the device and turn on the power. Use the original power cable.
- USB port: for connecting to LIS.
- Port for connecting an automatic pipette: The system will automatically start the measurement by pressing the pipette button.
- RS232 port: connection to a PC for data transfer.

Pic. 3-2 Configuration of the back panel of the analyzer "Technology Solution 4»

Caution: use the original power cord, do not use other power cables. Do not restart the device if you wait less than 10 seconds after turning off the device to avoid damage.

3.2 Specifications

- Measurement method: double magnetic circuit and magnetic ball.
- Four channels can run different tests simultaneously.
- Color touch-sensitive LCD display provides easy operation.
- Combinations of specific user tests.
- Intelligent temperature control, automatic tests, automatic timer and automatic mixing.

- Maintaining data quality control of PT, APTT, TT, FIB for 12 months.
- Open calibration curves for coagulation tests.
- Adding reagents to the cuvettes from the kits.
- Print results using the built-in printer.
- Data storage of 4000 samples.
- The software allows you to send analysis reports to the control computer.
- Connecting to LIS.

4. Functions

4.1 Functional scheme of the device

Tests

Choice of the test

- Test
- • Change of ID of the patient
- Printing of results

Management of reports of analyses

- Printing of results
- • Delete the selected report
- Delete all reports
- To unload reports

Parameters

- Test selection
- Information about the test

System settings

- Date and time
- Language setting
- Printing setting

Calibrations

- Test selection
- Changing the calibration curve

Quality control

- Test selection
- Diagram of quality control

Technical maintenance

- Password entry

Technical support

- Information about the device

4.2 Function buttons

	Moving the cursor up		Moving the cursor down
	The movement of the cursor to the right		The movement of the cursor to the left
	Page up		Page down
	Save and return		Edit the selected area
	Add information		Return back
	Printing an empty piece of paper		Printing data
	Return to test selection		Confirm the current operation and return to the previous menu
	Delete the current item		Return to the previous interface
	Selecting multiple tests to run simultaneously		Move the cursor one character to the left
	Timer in standby mode (digit color black)		Timer in the countdown process (the color of the digits is red)
	Timer – the countdown is completed		

5. Operation of the device

5.1 Installation

- Check all components according to packing list
- Remove all packaging materials
- Install the device on a flat table.

Caution: Roughness of the table may affect the accuracy of the test results..

- Attach the automatic pipette to the analyzer.
- Connect the device to the power source using the power cord.

Warning: Use a grounded socket.

If the difference between the temperature of the device and room temperature is more than 10 °C, please do not turn on the device immediately. It should be at room temperature for 24 hours before turning on.

5.2 Turning on

Turn the power button to the "I" position, the logo will appear on the LCD, and a beep will sound. If the printer is "enabled" by default, it will print the message "Technology Solution 4 semi-automatic coagulometer Report".

5.3 Beginning of use

Turn the power button to the "I" position and the LCD will display the logo (Pic.5-1), a beep will sound. If the printer is "enabled" by default, it will print the message "Technology Solution 4 semi-automatic coagulometer Report". Then the device automatically heats up. The warm-up time depends on the room temperature.

Pic. 5-1 Logo

5.4 Get zero point values


Tap anywhere on the screen and the system will automatically collect zero information if no cuvettes are insert in the measuring channels. The LEDs of the four measuring

channels will turn on together once, and a beep will sound.

Attention: please remove all cuvettes from the positions of the measuring channels and make sure that there are no foreign objects.

5.5 Choice of the test


After receiving the zero point values, the test selection menu appears (Pic. 5-2). In this menu, the user can select a test to test or set tests to run simultaneously.

Pic. 5-2 Choice of the test
— Select one or more tests to run.

5.5.1 Choice of one test


Select a test on the monitor in accordance with Pic. 5-2 and press "Enter", you will enter the analysis execution menu (Pic. 5-3).

Pic. 5-3 One Test

5.5.2 Choice of several tests

Select the "Multiple" item on the monitor in accordance with Pic. 5-2, you will enter the test group selection menu (Pic. 5-4). User can select or edit a group of frequently used tests.


Pic. 5-4 Group of selected tests
(CH1-CH4 - names of measurement channels)

5.5.2.1 Editing group of tests

Click "Edit" on the monitor in accordance with Pic. 5-4 to enter the editing menu (Pic. 5-5). You can add a test, modify or delete the entire test group.

Select the position that you want to change in Pic. 5-5 (the selected position shows a dark blue color). Click "Add" to enter the "Add Test" menu (Pic. 5-6). Click "Delete" to remove this test.

Select the tests that you want to add in Pic. 5-6 (the selected item shows a dark blue color), press "Enter" to save the changes.

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
Pic. 5-6 Adding Tests

When editing the test group is complete, click Back to return to the test group selection interface.

5.5.2.2 Choice of group of tests

Select a test group from the Multiple Tests menu. In Pic. 5-4, press "Enter", you will enter the measurement menu of several tests (Pic. 5-7). If the group's tests are not completed, channels to which the test is not define will show the status is prohibited.


Pic. 5-7 Measurement of multiple tests

5.6 Preparation for measurement

5.6.1 Changing the ID number

On the measurement menu, click the ID number area on different channels, you can enter the ID number change menu, the changed area will be highlight in red (for example, a separate test in Pic. 5-8). The user can manually change the identification number using the buttons.

When performing a test of one element, the identification number of each channel is independent, changing one ID number does not affect the other channel identification numbers. But when performing a test with several tests, the identification number is the same for all channels, one change in ID will lead to a change in all channel identification numbers.


Pic. 5-8 Changing the ID number


5.7 Preliminary incubation of samples

5.7.1 Beginning of a preliminary incubation

Click the clock icon below to start pre-incubation for different channels (Pic. 5-9). Four timers are design for four pre-incubation positions. When the pre-incubation begins, the time will be count down, when the pre-incubation is over, an audible signal will sound. To complete the incubation, click on the timer icon (Pic. 5-10).

Pic. 5-9 Pre-incubation


Pic. 5-10 Completion of pre-incubation

5.8 Test

The method for measuring multiple tests and with one test is similar. In this chapter, the measurement of one test is given as an example to explain the measurement procedure.

5.8.1 The beginning of the measurement

Using an automatic pipette, promptly add reagents to the test cells that are located in the measuring cells, the test will start for the corresponding channel and an audible signal will sound. The pre-incubation timer for the tested channels will disappear, the "Scroll" icon will disappear, and the "Test" icon will be replaced with the "Stop" icon (Pic. 5-11). Click "Stop" to cancel the current test, and click "Back" to return to the initial test screen.

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
5.8.2 End of the test


5.8.2.1 Normal end of the test

When the test is completed, "Completed" will appear on the screen and the test time will be shown "T = ..." and other test indicators (Pic. 5-12), the system will automatically save the test results. You can transfer the results to your PC or LIS system via the USB port. If the printer is ready, you can click the Print button to print the test results, and then return to the test waiting menu. If you do not need to print, click the Back button, the test results will not be printed, and you will return to the test waiting interface. Click the "Test" button, the test results will not be print, and you will return to the "test Selection" menu (Pic. 5-2).

Pic. 5-12 Completed PT Test


If the "FIB" test is being measured, the system will display the concentration value "C" = ... "after the time "T = ...". see Pic. 5-13.

Pic. 5-13 Completed test (FIB)

5.8.2.2 Abnormal end of the test

If there are no cuvettes in the measuring channels or there are other reasons that do not allow correct measurement, the screen for this channel will display "Completed", "No value". There will be no measurement result (Pic. 5-14), and the system will not save the test results.


Pic. 5-14 Completed test

Note: results can be printed after all tests are completed.

5.9 Main menu


The main menu is shown in Pic. 5-15, section 4.1 "Functional scheme of the device" of the function of each icon.

Pic. 5-15 Main menu

5.10 Management of reports


Click on the "Reports" icon to enter the report management menu. Select the appropriate option and click on the icon to enter (Pic. 5-16).


Pic. 5-16 Management of reports

5.10.1 Printing of reports


Click "Print reports" to enter the print menu (Pic. 5-17). In this menu, find the ID number that you want to print using the arrow buttons, and click the "print" button to search for it (Pic. 5-18). If the result for the ID number does not exist, the screen will display "no data found in the system". (Pic. 5-19). If results for an ID number are found, all tests for that ID number will be printed.


Pic. 5-17 Printing of reports


Pic. 5-18 Search of reports

Pic. 5-19 Can't find the report


Pic. 5-20 Printing of the report

Pic. 5-21 Completion of printing of the report

5.10.2 Deleting a single report

Click "Delete report" to enter the menu. Find the ID number that you want to delete using the arrow buttons, click Enter to delete it.


Pic. 5-22 The deletion is performed

After successfully deleting the saved data, the screen displays "Deleting reports completed", as shown in Pic. 5-22.


If there is no data for the ID number, the message "no data found in the system" will appear on the screen.

5.10.3 Deleting all reports


Click the "Delete all reports" button to enter the menu (Pic. 5-23). Click Enter to delete all reports. Deletion is complete when the progress bar shows 100 % (Pic. 5-24, 5-25).

Pic. 5-23 Deleting all reports

Pic. 5-24 Delete progress


Pic. 5-25 The deletion is performed

Note: do not click the "Return" button when deleting, otherwise it may cause a data error.

5.10.4 Uploading reports

Change the ID number by clicking the arrow buttons and upload the data for the desired ID number. Click Enter to upload all data for this ID number. When the upload process is complete, a message will appear on the screen, as shown in Pic. 5-26.


Pic. 26 Uploading data

5.11 Parameters

5.11.1 Choice of parameter

Click "Options" in the main menu to enter the menu (Pic. 5-27).

Pic. 5-27 Parameters

5.11.2 Setting of parameters

In the "Parameters" section, the user can set parameters. The parameters for the most used four tests (for example) as shown in Pic. 5-31 are shown below.

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
: moving the cursor to the left;

49

: use the "Up" and "Down" buttons to enter parameter data from 0 to 9 in the current row; when the cursor remains in the "Unit" parameter, click the button to select the desired block. For the FIB test, the value of "Upper limit of norm" and "lower limit of norm" will change accordingly when you select different units of measurement. For example, if the units are "g/l", "the Upper limit of normal" will be 4.00 and "low-normal" 2.00; if the unit "mg/DL", "Upper limit of normal" will automatically change to 400, and "low-normal" 200.  : Save to complete the settings. The APTT and TT parameters are return to the "Setting parameters" menu after setting.

5.12 System settings


Click "System settings" to enter the menu (Pic. 5-32).

Pic. 5-32 System settings

5.12.1 Date and time

Click the "Date and Time" button to set the date and time (Pic. 5-33).


Pic. 5-33 Setting the date and time

Click the arrow buttons to set the date and time, and click Save to save your changes.

5.12.2 Choice of language

Click "Select language" to select the interface language (Pic. 5-34). Click Save to save your changes.

Pic. 5-34 Choice of language

5.12.3 Print settings

Click "Print settings" to enter the menu (Pic. 5-35).

"Printer switch ""Enable" means that printing is enabled, "turn Off" means that the printer is turned off. In the "PT Options" section, you can choose which parameters to print when printing the PT report. Click Save to confirm the change.


Indicators on the screen display:

ПВ - the prothrombin clotting time;

МНО – the international normalized relation;

по Квику - prothrombin index for Quick;

ПТИ – prothrombin index.

Pic. 5-35 Print setting

When you print out the results of the prothrombin test (PT), the indicators will be displayed in the following form (Pic. 5-36):

```
* * * * *
Date: 2019-12-17
Item: PT
Ref.: 10.0--18.0    Sec
* * * * *

Sample      Test
201912170033
          T = 13.7 Sec
          INR = 1.11
          PTA = 72.8 %
          PTR = 1.11
```

Pic. 5-36 Example printouts of the results of the prothrombin test (PT)

The prothrombin test (PT) indicators are display on the printout

T - the prothrombin clotting time (PT when setting parameters);


INR – international normalized ratio (INR when setting parameters);

PTA - prothrombin index for Quick (по Квику - when setting parameters);

PTR – prothrombin index (ПТИ - when setting parameters).

5.13 Calibration

To perform calibration, enter the "Calibration" menu item and select the appropriate procedure (the selected item shows dark blue color), press "Enter".


Pic. 5-37 Entering the calibration menu

5.13.1 Set the calibration curve for the FIB test

The calibration values for the FIB curve are enter as follows:

The values from the passport to the set of calibrators are enter in the concentration column of calibrators. In the clotting time column, the values obtained when determining this indicator for the corresponding calibrator are enter. See Pic. 5-38.

Pic. 5-38 Values of 5 calibration points of the FIB test (units rev .: g / l)

: Click save to enter the FIB curve display menu (Pic. 5-38).

Click to return to the Calibration menu.

Function buttons


: moving the cursor to the right;

: moving the cursor to the left;

: Use the Move Up and Move Down buttons to enter parameter data from 0 to 9 in the current line.

Pic. 5-39 FIB test calibration curve

On the FIB curve menu, the x axis is concentration, the y axis is time, and the correlation coefficient "r" is shown in the middle. The closer the "r" value is to -1, the better the ratio of 5 points and vice versa.

Press to print the values of 5 calibration points (the printer must be turned on) (see Pic. 5-40).

Press  to return to the previous menu.


```
*****
Date:2019-08-06
Time:01:52:00
Item Name:FIB
Unit:g/L
C1: 6.56   T1: 7.5
C2: 3.28   T2: 14.0
C3: 1.64   T3: 28.0
C4: 1.09   T4: 37.8
C5: .82    T5: 53.1
a=-0.9318
b= 3.497
r=-0.9991
```

Pic. 5-40 Printed values of FIB calibration curve points

5.13.2 Set the calibration curve for the PT test

The calibration values for the PT curve are entered as follows:

Values from the passport to the plasma calibrator set and subsequent dilutions in percent by Quick are enter in the concentration column of calibrators. In the clotting time column, the values obtained when determining this indicator for the corresponding calibrator are enter. See Pic. 5-41.

Pic. 5-41 Values of 4 calibration points of the PT test

: Click save to enter the PT curve display menu (Pic. 5-41).

Click to return to the Calibration menu.

Function buttons


: moving the cursor to the right;

: moving the cursor to the left;

: Use the Move Up and Move Down buttons to enter parameter data from 0 to 9 in the current line.

Pic. 5-42 Calibration curve of the PT test

On the PT curve menu, the x axis is % Quick prothrombin according to Quick, the y axis is time, and the correlation coefficient "r" is shown in the middle. The closer the "r" value is to -1, the better the ratio of 4 points and vice versa.

Click  to print the calibration points (the printer must be turn on). See Figure 5-40.

Press  to return to the previous menu.

```

*****
Date:2019-08-06
Time:01:53:00
Item Name:PT
Unit:Sec
C1: 100      T1: 12.2
C2:  50      T2: 17.3
C3:  25      T3: 24.8
C4:  13      T4: 40.4
a=-0.5792
b= 2.230
r=-0.9952

```

Pic. 5-43 Printed PT calibration curve points

5.14 Quality control

Click "Quality Control" in the main menu to enter the menu (Pic. 5-44).

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
→ : moving the cursor to the right;

 : moving the cursor to the left;

  : Use the Move Up and Move Down buttons to enter parameter data from 0 to 9 in the current line.

5.14.1 Representation of quality control on a graph

Click  to enter the "View Quality Control on the Graph" menu (Figure 5-44). The QC chart of the current month by default will be display on the screen. Each test can store a maximum of 12-month QC data. If the test value is "0", it does not display a point. In Figure 5-44, the X axis is the date axis (unit: day), and the Y axis is the test time axis (unit: seconds). The middle line of the X axis is the target value, the two horizontal lines of the upper and lower level of the $x + 2SD$ and $-2SD$ axis, the upper and lower horizontal lines of the $X + 3SD$ and $-3SD$ axis. When the QC is greater than $+2S$ or less than $-2S$, the user must check the instrument, reagent, or methodology to analyze the reasons for the deviation; when the QC is greater than $+3S$ or less than $-3S$, the user must stop using this instrument and analyze the reasons for the deviation.

Pic. 5-45 Representation of quality control on a graph

Function buttons

: Presentation of QC graph of last month. If the interface is the current month, the system will display the QC value 12 months ago.

: Next month QC graph view, 12-month cycle


: To return to the QC menu.

5.15 Service menu

Login is done using a password, and access is provide by a qualified engineer.

5.16 Technical support

Select "Technical Support" in the main menu, the analyzer's name, analyzer serial number, software version, service hotline, website and the name of the company that supplies the devices in the Russian Federation will appear on the screen (Pic. 5-46). Touch any point on the screen to return to the main menu.

Pic. 5-46 Technical support

5.17 Printer

The device uses a built-in thermal printer. Thermal paper 58 mm wide is required.

Attention. The printer cannot print on normal printing paper.

5.17.1 Installing paper in the printer

The printer will not print if there is no paper to print or if the printer cover is not closed well enough. You must install the print paper manually. Press the button to open the printer cover, place the print paper in the correct direction in the printer head, and then close the printer cover tightly.

5.17.2 Removing paper from the printer

Press the button to open the printer cover, then you can remove the printing paper.

5.17.3 Thermal printer maintenance

- Use high-quality thermal paper; paper quality affects the clarity of printed characters.
- Regularly clean the thermal roller from paper dust.
- Use a cotton swab to gently clean the heating pad of the thermal printer.

5.18 Connection with PC

To connect the device to a PC for data transfer, use the "USB port" on the back of the device.

6. Work Precautions

- 1) The power supply voltage of the device is 220 V. If the network parameters do not allow to fulfill this requirement, it is necessary to use an uninterruptible power supply of alternating current with a power of more than 200 watts.
- 2) Do not use the analyzer in a room containing flammable substances to avoid explosion.
- 3) Installation and use of the device should be carry out away from a strong electromagnetic field. When working with the device, do not use high-frequency equipment (mobile phone, etc.). High-frequency radiation can cause the device to malfunction and affect the measurement results.
- 4) The analyzer must be operate by trained personnel. The supplier company during installation provides training in the operation of the device.
- 5) **The user must strictly follow the operating instructions for this analyzer.**
- 6) To ensure the high quality of the results obtained, only those reagents and test cuvettes recommended by the organization supplying the instruments in the Russian Federation should be use.
- 7) **The user must strictly follow the instructions attached to the reagents.**
- 8) Please do not spill test samples to prevent biological contamination.
- 9) Remove all samples and reagents when moving or repairing the analyzer to avoid contamination or potential biological contamination.

This analyzer was thoroughly tested and configured before leaving the factory.

Do not change the settings of parameters and accessories yourself.

7. Important notes and maintenance

7.1 Important notes

- 1) The preheating of the analyzer is 30 minutes (heating to 37.0 °C, at an ambient temperature of 20.0 °C). Please do not operate the instrument until the temperature is stable, to ensure accurate measurement.
- 2) When placing the device, avoid direct sunlight and high brightness of lighting, strong heat and strong electromagnetic interference.
- 3) Please do not touch the bottom half of the test cuvette to ensure accurate measurement.
- 4) The resource of an automatic pipette is at least 5000 times of use. After this number of uses, it is recommend to replace it.
- 5) Please control the force by pressing the touch screen, we also recommend using the touch pen.
- 6) The print head is easily corroded, please keep it dry and clean.
- 7) Remove the test cuvettes from the measuring cells immediately after the analysis and keep the measuring cells clean after the tests are completed.
- 8) Please do not touch the screen during the test to avoid interruption of the test and data loss.
- 9) The temperature information will not be update during the test.
- 10) It is recommend to regularly monitor the system time on the device due to the possibility of fluctuations in the system clock.
- 11) Pay attention to the series and expiration date of reagents. If a set series change occurs, calibrate the procedure (for calibrated tests).
- 12) It is proposed to use those consumables to ensure the accuracy of the results of clinical trials that are recommended by the organization that supplies the devices in the Russian Federation
- 13) Specialized disposable blood plasma cuvettes and reagents formed during testing of the sample should be processed and disposed of in accordance with the norms established by law, their reuse and uncontrolled disposal are prohibited.

Waste from the production of the device belongs to class 5 according to the degree of danger to the environment. Waste is not harmful to the environment. Measuring ditches, not used for their intended purpose, belong to class A -

epidemiologically safe waste, similar in composition to municipal solid waste according to SanPiN 2.1.7.2790-10.

Ditches measuring after use in medic. institutions belong to class B - epidemiologically hazardous waste according to SanPiN 2.1.7.2790-10.

7.2 Maintenance

- 1) Keep the appliance clean. Use a soft cloth to remove dust from the surface of the analyzer. To clean a medical device, an aqueous solution containing a neutral detergent with the following characteristics is used: less than 5% anionic surfactants, less than 5% nonionic surfactants, less than 5% amphoteric surfactants (for example, CIF "Active Lemon", manufacturer of Unilever Rus LLC or similar), to remove contamination from the surface of the analyzer. Detergents containing organic solvents and abrasives are not used.
- 2) This device is a high-precision measuring device, it should be installed on a flat surface, after use cover it with a dust-proof cover.
- 3) To ensure accurate test results, calibrate the analyzer regularly.
- 4) Routine and service maintenance is not provided.

8. System errors and their elimination

Errors	Check of the reasons	Decision
After starting the display does not work	Check the power supply. Check the fuse.	Reconnect the power cord or replace the fuse.
The printer does not print.	The printer is set to "off" or the thermal paper has run out.	Set the printer to "on" or insert paper for printing.
The color of the printed characters is light.	Fault of thermoprinting.	Replace the thermal paper.
The printer prints without paper.	The print head is dirty.	Use a brush or dry cloth to clean the dust on the printhead.
The device does not respond to a tap on the screen	The click position is incorrect. The loss of information about position of the screen. The device tests samples. The device saves parameters. The printer prints.	Tap the correct place on the display. Calibrate. Wait until the analyzer finishes working.
The test result is not accurate, or the spread of indicators is too large.	There are bubbles in the test cuvette. The test cell is dirty.	Check if there are any bubbles on the inner wall of the cell. Replace the cuvette.

9. Manufacturer's warranty

The lifetime of the device is 5 years

The authorized representative in the Russian Federation of LLC "Laboratory Solutions" provides the warranty life of the device 3 years from the date of its commissioning, provide that the operating rules given in this manual are observed.

The operator shall use the instrument in accordance with this operating manual.

For correct use of the device during commissioning, the representative of the authorized representative in the Russian Federation of LLC "Laboratory Solutions" will conduct training on operation on the device.

If there is a fault in the instrument, contact to qualified specialists of "Laboratory Solutions" LLC.

In order to maintain the warranty for the device, it is necessary to use reagent kits and cuvettes recommended by the authorized representative in the Russian Federation of LLC "Laboratory Solutions" during the warranty period.

In addition, we pay attention to the user of the device that the possibility of operating a medical product of one manufacturer together with consumables (sets of reagents, cuvettes, one-time ball) of another manufacturer is determined by an authorized representative in the Russian Federation of LLC "Laboratory Solutions."

In case of violation of the above conditions, the authorized representative in the Russian Federation of LLC "Laboratory Solutions" removes any responsibility for the quality of the obtained results and for the operability of the device, and also has the right to refuse replacement on guarantee of the failed device or any part thereof.

Routine and service maintenance is not provide. In case of failure of any unit of the device, contact specialists who have a permit certificate for this model of the device.

The permit certificate is issue by an authorized representative in the Russian Federation of LLC "Laboratory Solutions."

In case of failure due to incorrect operation of the device, as well as failure to comply with the above mention conditions, the warranty is terminate.

10. Packaging, storage and transportation

The device is pack in a box of corrugated cardboard, ensuring the safety of the medical device from climatic and mechanical influences during transportation.

Inside the box, the device is surround on all sides by foam gaskets.

It is necessary to store the product in dry rooms at a temperature of $0^{\circ}\text{C} + 40^{\circ}\text{C}$ and relative humidity $\leq 90\%$

Keep away from moisture.

Avoid heating or storage in direct sunlight, dust or aggressive gases.

The research room should have an operating temperature of $+10^{\circ}\text{C}$ to 30°C and relative humidity $\leq 70\%$

If the appliance cools down during transport, please wait while it warms up to room temperature before turning it on.

When storing, the product must be pack in the packaging container of the supplier (manufacturer).

It is allowed to transport all types of vehicles in covered vehicles in accordance with the rules of transportation applicable to this type of transport at a temperature of -30°C to $+55^{\circ}\text{C}$ and relative humidity up to 90%.

The device should be transport only in the packaging of the supplier (manufacturer).

It is allow to transport the product in additional transport packaging of the transport company.

The method of laying transport containers with products should exclude the possibility of their movement.

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