Tre Hargett  
Secretary of State

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Tre Hargett  
Secretary of State

**Division of Business Services  
Department of State**

State of Tennessee  
312 Rosa L. Parks AVE, 6th FL  
Nashville, TN 37243-1102

**APOSTILLE**

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1. Country: United States of America
2. This public document has been signed by SHERRY WITT
3. acting in the capacity of COUNTY CLERK OF KNOX COUNTY, TENNESSEE
4. bears the seal/stamp of COUNTY SEAL OF KNOX COUNTY, TENNESSEE

**Certified**

5. at Nashville, Tennessee
6. the 11th day of February, 2020
7. by TRE HARGETT, SECRETARY OF STATE, STATE OF TENNESSEE
8. No. 20-02156

**10. Signature**

*Tre Hargett*  
Secretary of State

by *Doreen Hackney*

To verify authenticity, visit <http://tnbear.tn.gov/apostille/verify.aspx> and enter the document # (line 8 above) from this certificate. For further information please contact the Business Services Division at (615) 741-0536.

**SHERRY WITT**  
KNOX COUNTY CLERK  
300 MAIN STREET  
KNOXVILLE, TN 37902  
865-215-2385

## Notary Public Authentication

I, SHERRY WITT, County Clerk of the County of KNOX, Tennessee do hereby certify that ANNA M KIRK was duly elected and commissioned a Notary Public At-Large for the State of Tennessee on 2018-11-27, and is recorded in Book 109 Page 163 in this office. The above-mentioned Notary Public At-Large commission does not expire until 11/27/2022.

I further certify that to the best of my knowledge, information and beliefs, the signature affixed to the attached document is the true and genuine signature of ANNA M KIRK, Notary Public for the State of Tennessee.

In Witness Whereof, I have here to affixed my Signature and the seal of the County of KNOX this day.

Sherry Witt  
KNOX COUNTY CLERK

DAVID GARRETT  
Deputy Clerk

David Garrett

February 07, 2020  
Date



This is to confirm that this document has been  
printed by Siemens Medical Solutions USA, Inc.,  
2501 N. Barrington Rd, Hoffman Estates, IL 60192 USA.  
Tabitha Estes

## User Manual addendum for Russia

*System of combined single-photon emission and X-ray computed tomography of  
Symbia T6 with accessories*

On this 7th day of February 2020, I attest that the attached document is a true, exact, complete, and unaltered copy of the Artwork of the User Manual addendum for Russia

  
Tabitha Estes – Regulatory Affairs Specialist

  
Notary Public

A Notary Public in and for the County of Knox, State of TN.

## 1. The name of the Medical Device.

System of combined single-photon emission and X-ray computed tomography of Symbia T6 with accessories.

## 2. Manufacturing information of the Medical Device.

Legal manufacturer:

Siemens Medical Solutions USA, Inc.  
Molecular Imaging  
2501 N. Barrington Road,  
Hoffman Estates, Illinois 60192  
USA  
Tel.: +1 847 304 7700  
Email: [www.siemens.com/healthcare](http://www.siemens.com/healthcare)

Manufacturing sites:

Siemens Medical Solutions USA, Inc.  
Molecular Imaging  
2501 N. Barrington Road,  
Hoffman Estates, Illinois 60192  
USA  
Tel.: +1 847 304 7700  
Email: [www.siemens.com/healthcare](http://www.siemens.com/healthcare)

The name and address of the local representative in Russia:

Siemens Healthcare LLC  
Dubinskaya street 96  
115093 Moscow,  
Russia  
Tel.: +7 495 737 10 00  
Fax: +7 495 737 13 20  
[www.healthcare.siemens.ru](http://www.healthcare.siemens.ru)

### 3. Designation of the Medical Device.

#### 3.1. Indications

Symbia T6 combined single photon emission and X-ray computed tomography system with accessories is intended for use by appropriately trained health care professionals to aid in detecting, localizing, diagnosing, staging and restaging of lesions, tumors, disease and organ function for the evaluation of diseases and disorders such as, but not limited to, cardiovascular disease, neurological disorders and cancer. The images produced by the system can also be used by the physician to aid in radiotherapy treatment planning and interventional radiology procedures.

**SPECT:** The SPECT component is intended to detect or image the distribution of radionuclides in the body or organ (physiology), using the following techniques: planar imaging, whole body imaging, and tomographic imaging for isotopes with energies up to 588 keV.

**CT:** The CT component is intended to produce cross-sectional images of the body by computer reconstruction of x-ray transmission data (anatomy) from either the same axial plane taken at different angles or spiral planes taken at different angles.

**SPECT+CT:** The SPECT and CT components used together acquire SPECT/CT images. The SPECT images can be corrected for attenuation with the CT images, and can be combined (image registration) to merge the patient's physiological (SPECT) and anatomical (CT) images.

**Software:** The syngo MI Applications software is a display and analysis package intended to aid the clinician in the assessment and quantification of pathologies in images produced from SPECT, PET, CT, and other imaging modalities.

#### 3.2. Intended for use

The Symbia T6 is radiological imaging system that combine a single photon emission computed tomography (SPECT) camera system for nuclear medicine images, and a computed tomography (CT) camera system for x-ray images.

The SPECT system is intended to detect the location and distribution of gamma ray radionuclides in the body and produce cross-sectional images through computer reconstruction of the data, and the CT system is intended to produce cross-sectional images of the body by computer reconstruction of x-ray transmission data from the same axial plane taken at different angles. The nuclear medicine images and the x-ray images may be registered and displayed in a fused format (overlaid in the same orientation) for the anatomical localization of the nuclear medicine data (that is, distribution of radiopharmaceuticals).

The SPECT and CT portions of the system may be used independently or in combination, and may include signal analysis and display equipment, patient and equipment supports, radionuclide anatomical markers, component parts, and accessories. The SPECT and CT images may be transferred to other systems for radiation therapy planning or additional uses.

### **3.3. Contraindications:**

There are no known specific situations that contraindicate the use of the system.

### **3.4. Possible side effects:**

In clinical use, CT systems produce X-rays which are a type of ionizing radiation. Ionizing radiation has the potential to cause long-term biological effects in the human body such as cancer, fetal effects, and genetic effects.

### **3.5. Conditions for use:**

The Symbia T6 system must be operated according to the indication and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training.

### **3.6. Staff qualifications**

Only properly trained and authorized users with specialized knowledge such as trained physicians, trained radiologists, or trained technologists should operate the system and its accessories. Improper or unauthorized use of the system can cause injury to the user or patient, damage to the system, or incorrect diagnosis or treatment due to misinterpretation of image information. Nuclear medicine training, medical training, and specific training on the Symbia T6 are mandatory.

### **3.7. Special treatment consideration**

Certain atypical situations require special consideration by the medical professional on behalf of the patient. These situations include pregnancy, breastfeeding, pediatric cases and adults with chronic diseases. In these cases the medical professionals discuss with each patient the benefits of research while informing them that the risks are higher than in an average patient, and there are standard techniques for minimizing risks for these patients. This is also true in the case of patients with implanted medical devices.

### **3.8. Operation conditions (humidity, temperature, pressure)**

| Parameter           | Range value                  |
|---------------------|------------------------------|
| Temperature         | 66°F (+18°C) to 86°F (+30°C) |
| Relative humidity   | 20% to 80 % non - condensing |
| Barometric pressure | 700 -1060 hPa                |

## 4. Classification of the Medical Device.

### 4.1. Potential risk class:

Classification according to Annex IX, MDD – Class IIb

### 4.2. Type of the Medical Device in accordance with the Russian Nomenclature classification:

280730

### 4.3. Class of the Medical Device in accordance with IEC 60601-1:

Class I according to IEC 60601-1:2005;

Type B applied part (pallet)

Defibrillator Proof Type BF applied part (ECG trigger)

### 4.4. Type of contact with the human body:

The Medical Device is non-sterile, nonexpendable.

Device surfaces that can come into contact with patients are evaluated and proven to be biocompatible according to ISO 10993:2011.

### 4.5. Software Risk classification / Level of concern:

Risk Classification according to IEC 62304 – Class B.

Level of concern according to FDA Software Guidance – Moderate Level of Concern.

Version SW – VB20

Name of SW -- syngo® MI Applications.

#### **4.6. Laser safety classification**

Laser safety classification acc. to EN 60825-1 – Class 1M.

#### **4.7. Classification according to the degree of protection against ingress of water as detailed in the current edition of IEC 529**

Classification according to the degree of protection against ingress of water as detailed in the current edition of IEC 529 – IPX0 (ordinary).

#### **4.8. Classification according to the mode of operation**

Continuous operation with intermittent loading.

#### **4.9. Classification according to SanPiN 2.1.7.2790-10**

Class A.

### **5. Operation of the Medical Device.**


#### **5.1. Physical principles:**

Single photon emission computed tomography (SPECT) is a type of emission tomography; diagnostic method for creating tomographic images of the distribution of radionuclides. In SPECT, radiopharmaceuticals labeled with radioisotopes are used, the nuclei of which emit only one gamma quantum for each act of radioactive decay.

SPECT / CT (single photon emission computed tomography combined with computed tomography) is a highly sensitive diagnostic method that allows you to scan the entire body, a limited area of the body or a separate organ.

The image of metabolic disturbances obtained on SPECT is superimposed on computed tomography (CT), which makes it possible to see a three-dimensional reconstruction of the processes and anatomical changes. Thus, with SPECT / CT examination, it is possible to obtain information about a human disease both at the molecular and cellular level, and at the organ level.


## 6. The overview of the Medical Device.

Main parts of Symbia T combined single photon emission and X-ray computed tomography system with accessories:


|    |                                               |
|----|-----------------------------------------------|
| 1  | Emergency Stop Button                         |
| 2  | Hand Control                                  |
| 3  | Nuclear Medicine (NM) Subsystem               |
| 4  | Computed Tomography (CT) Subsystem            |
| 5  | Rear Bed                                      |
| 6  | Siemens Nuclear Acquisition Controller (SNAC) |
| 7  | Power Switches                                |
| 8  | Patient Positioning Monitor (PPM)             |
| 9  | Touch Pads                                    |
| 10 | NM Detectors                                  |
| 11 | Integrated Source Holder (retracted)          |
| 12 | Pallet                                        |

|                                                                                     |                                                  |
|-------------------------------------------------------------------------------------|--------------------------------------------------|
| 13                                                                                  | Patient Bed                                      |
|  |                                                  |
| Main parts of patient bed                                                           |                                                  |
| 1                                                                                   | Rear Bed                                         |
| 2                                                                                   | Foot End                                         |
| 3                                                                                   | Pallet with Pad                                  |
| 4                                                                                   | Patient Bed Handles                              |
| 5                                                                                   | Docking Pin (beneath bed) - Foot Pedal not shown |
| 6                                                                                   | ECG Connection                                   |
| 7                                                                                   | Pallet Handle                                    |
| 8                                                                                   | Bed Index                                        |
| 9                                                                                   | Bed Index Handle                                 |
| 10                                                                                  | Head End                                         |


## 7. List of components and accessories of the Medical Device:

| Название принадлежности             | Name of Component           | Description                                                                                                                                                                                                                                                                                                                                                                                | Photograph                                                                          |
|-------------------------------------|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| I. Блок базовый Symbia T6:          | Symbia T6 (basic block):    |                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                     |
| 1. Гентри ОФЭКТ/КТ интегрированное. | Integrated SPECT•CT Gantry: | The gantry has two Variable Angle SPECT detectors and an open design with a 70 cm (27.6 in.) patient opening. The unobstructed gantry base permits planar imaging of seated and standing patients and patients on wheelchairs, or on standard imaging tables, gurneys and hospital beds. The Ultra Fast Ceramic multislice spiral CT detector rotates at 120 RPM (0.5 sec per revolution). |  |

|                                               |                 |                                                                                                                                |                                                                                      |
|-----------------------------------------------|-----------------|--------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 2. Детектор высокоразрешающий 3/8, 5/8 дюйма. | SPECT Detectors | SPECT Detector detects photons                                                                                                 |   |
| 3. Коллиматор низких энергий многоцелевой.    | Collimators     | The collimator allows to define the direction of photon emission.                                                              |   |
| 4. Стол пациента.                             | Patient Table   | The patient bed supports the pallet and bed index, and raises them to go between the detectors, and through the gantry's ring. |  |

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
|                                                                                              |                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                                                                                       |
|----------------------------------------------------------------------------------------------|-----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <p>9. Модуль программно-аппаратный для модуляции тока на трубке в процессе сканирования.</p> | <p>CARE Dose 4D:</p>                                            | <p>Care Dose 4D automatically determines the minimal x-ray dose level needed to obtain optimal image quality, for all scan modes. The initial or starting tube current for every axial slice position is determined from the topogram image. Then, during the data acquisition for each axial slice, the x-ray attenuation values are closely monitored and the tube current is adjusted, on a real time basis, to optimize the x-ray dose level for the specific organs and anatomy in the x-ray path.</p> |    |
| <p>10. Монитор специальный цветной.</p>                                                      | <p>Monitor, 19" LCD<br/>DICOM<br/><br/>Monitor: 19 inch LCD</p> | <p>Color monitor</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |  |

|                                            |                            |                                                                                                                                                                                |                                                                                     |
|--------------------------------------------|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 11. Клавиатура, пульт управления КТ, мышь. | Keyboard<br>CT Control Box | The CT control box is used to start and stop a CT acquisition's radiation, move the patient bed, communicate with the patient, and move the patient bed's pallet.              |  |
| 12. Монитор позиционирования пациента.     | Monitor, PPM               | The patient positioning monitor (PPM) displays patient positioning, collimator changing, and detector reconfiguration information. It also displays an acquisition's progress. |  |

|                                  |                         |                                                                                                                                                                                            |                                                                                      |
|----------------------------------|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 13. Ручной пульт управления.     | Hand Controller         | The hand control is used to manually communicate with the system. It can control gantry, patient bed, and detector movements, as well as initiate acquisitions and programmed positioning. |   |
| 14. Опция удаленной диагностики. | Siemens Remote Services | A broadband VPN connection is required for full remote service functionality and optimal system uptime.                                                                                    |  |

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










|                                          |                                |                                                                                                                                                                                                                                                                                                                                 |                                                                                     |
|------------------------------------------|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 8. Тележка коллиматорная.                | Symbia Collimator Cart         | The collimator cart is designed to hold extra collimators and allows collimator exchange without removing the bed.                                                                                                                                                                                                              |  |
| 9. Опция удаленного контроля за сканером | syngo MI Remote Scanner Status | Remotely monitor the status of Symbia scanners at a single location within the department. Monitoring is supported by any device that has both access to the department's network and has web browser capability. The interactive status of an acquisition can be conveniently accessed from almost anywhere in the department. |  |

|                                                                                                |                                |                                                                                                                                                                              |                                                                                      |
|------------------------------------------------------------------------------------------------|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 10. Информационно-обучающая система e.media.                                                   | e.media                        | The e.media patient comfort and education package integrates high quality video and sound through the color touch screen patient positioning monitor.                        |   |
| 11. Лицензионные ключи для количественной оценки динамических КТ данных головного мозга.       | syngo Perfusion CT             | The syngo Perfusion CT option aids in the quantitative evaluation of dynamic CT data of the brain following the injection of a highly concentrated iodinated contrast bolus. |   |
| 12. Лицензионные ключи для мультимодальной визуализации и количественной оценки в кардиологии. | Cedars Fusion/Coronary Overlay | This option provides the ability to fuse a patient's coronary vessels onto the left ventricular SPECT or PET perfusion surfaces of the heart.                                |  |




|                                                                     |                   |                                                                                                                                                                                                                               |                                                                                       |
|---------------------------------------------------------------------|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 13. Лицензионные ключи для проведения онкологических исследований.  | Oncology Engine   | The Oncology Engine facilitates lesion detection by enabling the visualization, volumetric analysis, and fusion of SPECT studies as well registration of images from other, independently acquired modalities (e.g., CT, MR). |    |
| 14. Лицензионные ключи для количественной оценки перфузии.          | Cardiology Engine | The Cardiology Engine assists in the diagnosis and quantitative assessment of coronary artery disease by enabling the visualization of SPECT and/or CT studies as well as quantified perfusion assessment.                    |   |
| 15. Лицензионные ключи для проведения неврологических исследований. | Neurology Engine  | The Neurology Engine enables the quantification of neurology examinations.                                                                                                                                                    |  |



|                                                                                    |                        |                                                                                                                                                                                                                                |                                                                                       |
|------------------------------------------------------------------------------------|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 16. Лицензионные ключи для количественного анализа узелковых образований в легких. | CT Lung CAD            | syngo.CT Lung CAD is a computer-aided detection tool designed to assist radiologists in the detection of solid pulmonary nodules, partial-solid, and Ground-Glass Nodules (GGN) during review of CT examinations of the chest. |    |
| 17. Лицензионные ключи для виртуальной эндоскопии.                                 | Fly Through for Symbia | The syngo Fly-Through option provides high quality surface shaded (SSD) and volume rendering techniques (VRT) for virtual endoscopic viewing.                                                                                  |    |
| 18. Лицензионные ключи для проведения виртуальной колонографии.                    | CT Colonography        | syngo.CT Colonography combines the advantages of 2D and 3D reading strategies. Flexible screen layouts and dual monitor support permit instant switching between the 3D endoscopic                                             |  |

|                                                                              |             |                                                                                                                                                                                                                                                                                       |                                                                                     |
|------------------------------------------------------------------------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                                                              |             | view and the corresponding 2D images.                                                                                                                                                                                                                                                 |                                                                                     |
| 19. Лицензионные ключи для количественного определения КТ значений в легких. | CT Pulmo 3D | <p>CT Pulmo 3D is an application for clinical assessment of COPD (Chronic Obstructive Pulmonary Disease). Based on the user's preference, the application provides different automated 3D quantifications of a native chest CT scan for assessment of emphysema and lung airways.</p> |  |

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




|                                                                        |            |                                                                                                                                                              |                                                                                     |
|------------------------------------------------------------------------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                                                        |            | The IDL Runtime Software enables execution of custom applications. Custom applications can be developed and implemented by using the IDL Developer Software. |                                                                                     |
| 30. Лицензионные ключи для контроля поступления контрастного вещества. | CARE Bolus | The CARE Bolus option automatically triggers the start of a CT spiral scan when contrast media appears in sufficient density at a user defined location.     |  |

|                                                                         |                        |                                                                                                                                                                                                        |                                                                                       |
|-------------------------------------------------------------------------|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 31. Лицензионные ключи для обеспечения безопасности.                    | syngo Security Package | Software option providing enhanced security features including user management and audit trail functionality.                                                                                          |   |
| 32. Лицензионные ключи для объемной визуализации и анализа изображений. | FusedVision3D          | The FusedVision3D is a Volume Rendering Technique that provides visualization of fused anatomical and functional volumes via projection of the volumes onto an arbitrary oriented plane in full screen |  |


|                                                 |                  |                                                                                           |                                                                                     |
|-------------------------------------------------|------------------|-------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                                 |                  | mode or together with the 3-orthogonal fused datasets.                                    |                                                                                     |
| 33. Лицензионные ключи для обработки данных MI. | Organ Processing | Processing software package that provides cardiac and other organ-based SPECT processing. |  |


|                                                                         |                                |                                                                                                                                                                                                    |                                                                                       |
|-------------------------------------------------------------------------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 34. Лицензионные ключи для перевода изображений с ECAT в формате DICOM. | syngo ECAT to DICOM Translator | The ECAT to DICOM image translator converts native ECAT 7.2 image data to PET DICOM image data and automatically imports the DICOM images into the syngo MultiModality Workplace's syngo database. |   |
| 35. Стол рабочий для оператора.                                         | Computer Desk                  | The desk to accommodate the control components and color monitors.                                                                                                                                 |  |

|                                                     |                   |                                                                                                                                            |                                                                                      |
|-----------------------------------------------------|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 36. Стул для оператора с антистатическим покрытием. | Swivel Tilt Chair | Mid-pivot control mechanism reduces amount seat front rises during recline, to enhance circulation and reduce muscle fatigue to lower body |   |
| 37. Стекло рентгенозащитное (просвинцованное).      | Diagnostic Window | Diagnostic Window                                                                                                                          |  |


|                                   |                                     |                                                                                                                                                        |                                                                                      |
|-----------------------------------|-------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 38. Триггер ЭКГ.                  | ECG Gate w/out Strip Chart Recorder | ECG gating system provides display of ECG waveform and trigger indicator. Power and trigger signal connect conveniently the bottom of the Patient Bed. |   |
| 39. Кабель ЭКГ с электродами      | ECG cables                          | ECG cables                                                                                                                                             |   |
| 40. Электроцит распределительный. | Distributor                         | Distributor                                                                                                                                            |  |

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## 8. Additional Technical Parameters.

### 8.1. Minimum room size:

Assembly and Installation on site is performed by Siemens staff or Siemens trained and approved contractors.


Figure - Minimum Room Size

| Scanner room size   | 3.73 m (12 ft 3 in) x 6.32 m (20 ft 8 in) |                     |                                    |
|---------------------|-------------------------------------------|---------------------|------------------------------------|
| Ceiling height      | 2.44 m (8 ft 0 in)                        |                     |                                    |
| Hung ceiling height | 2.29 m (7 ft 6 in)                        |                     |                                    |
| Room Diagram Label  | Item Name                                 | Weight              | Heat Output                        |
| G                   | Symbia T Series gantry                    | 2369 kg (5,224 lb)  | 3413 BTU/h, 1.0 kW                 |
| CT                  | CT components                             | 1129 kg (2,490 lb)  | <1365 BTU/h, <0.4 kW <sup>16</sup> |
| T                   | Symbia T Series imaging table             | 950 kg (2,096 lb)   | –                                  |
| R                   | Symbia T Series rear PHS                  | 188.3 kg (415.2 lb) | –                                  |
| LCB                 | Line connection box                       | –                   | 1365 BTU/h, 0.4 kW                 |



| Control Room Heat output            |                                                                                                                                                                                                                                                             |
|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Acquisition computers               | 2389 BTU/h, 0.7 kW<br>(PMMA Phantom.<br>Absorbed dose for<br>reference material air.<br>Max. deviation $\pm 30\%$ .<br>Expected deviation $\pm 15\%$ .<br>Slice $>1\text{mm}$ . Please note<br>that these specifications<br>are CTDI <sub>100</sub> values) |
| Advanced reconstruction workstation | 2142 BTU/h, 0.8 kW                                                                                                                                                                                                                                          |

#### Workplace Proximity to Patient:

A perimeter distance of 1.83 meters (6.0 feet) must be maintained around the patient bed, with a distance of 2.5 meters (8.2 feet) between floor and ceiling.

## 8.2. Technical requirements. Key parameters & values provided by manufacturer. Characteristics (properties)

Maximum permissible deviation of technical characteristics no more than +/- 10%, unless otherwise specified

### Power Requirements

| SPECT input voltage | Single-Phase 200/208/220/230/240 VAC ~ 50/60 Hz                        |
|---------------------|------------------------------------------------------------------------|
| CT input voltage    | Three-Phase 380/400/420/440/460/480 VAC ~ 50/60 Hz                     |
| Electrical supply   | Short term input power 72.2 kVA<br>Continuous term input power 4.5 kVA |

### System Specifications

| System Software (Standard)                           |                                                |
|------------------------------------------------------|------------------------------------------------|
| Multimodality viewing software                       | Topogram, spiral, sequential acquisition modes |
| Planar (static) acquisition                          | Extended field of view (FoV)                   |
| Dynamic acquisition                                  | CARE Dose4D™                                   |
| Whole-body acquisition                               | SureView™                                      |
| Whole-body SPECT acquisition                         | Flexible kV (80, 110, 130 kV)                  |
| Gated acquisition                                    | DICOM structured dose report                   |
| Gated SPECT acquisition                              | syngo® archiving and network                   |
| Dynamic SPECT acquisition                            | Asynchronous reconstruction                    |
| Filtered back projection                             | Image display                                  |
| Flash 3D iterative reconstruction                    | Multiplanar reconstruction                     |
| Scatter correction                                   | Video capture and editing tool                 |
| CT attenuation correction                            | CT scan protocol assistant                     |
| Cardiac half-time imaging                            | WorkStream4D™ (direct 3D-reconstruction)       |
| Remote diagnostic services                           | syngo Viewing                                  |
| PPM display and interaction at acquisition workplace | syngo Filming                                  |
| Gated study beat normalization                       | syngo Dynamic Evaluation                       |
| Automatic and manual motion correction               |                                                |

| System Software (Optional)                                                  |                                        |
|-----------------------------------------------------------------------------|----------------------------------------|
| 3D measured collimator characterization for IQ•SPECT (hole, shape and size) | Neurology Engine<br>Scenium            |
| Organ Processing for Symbia                                                 | syngo media viewer                     |
| Planar half-time imaging                                                    | syngo CARE Bolus                       |
| syngo MI Remote Scanner Status                                              | syngo Perfusion CT                     |
| syngo Security Package                                                      | syngo Pulmo CT                         |
| Heartview CT                                                                | syngo Fly Through for Symbia           |
| Cedars Cardiology Engine                                                    | Multi-series CT attenuation correction |

#### SPECT Specifications

| Gantry Dimensions                                 |                                    |
|---------------------------------------------------|------------------------------------|
| Height                                            | 225 cm (7 ft 4.7 in)               |
| Width                                             | 231 cm (7 ft 7 in)                 |
| Depth                                             | 203 cm (6 ft 8 in)                 |
| Axis of rotation (from floor)                     | 104 cm (3 ft 5 in)                 |
| Weight                                            | 3,506 kg (7,714 lb)                |
| Min./max. patient opening (HE Coll)               | 12 cm (4.7 in)/65.4 cm (25.7 in)   |
| Min./max. patient opening (LEHR Coll)             | 19.2 cm (7.6 in)/72.6 cm (28.6 in) |
| Patient positioning monitor                       | 15" flat panel color LCD display   |
| Tunnel opening                                    | 70 cm aperture (27.6 in)           |
| Tunnel length                                     | 89 cm (35 in)                      |
| Distance between SPECT and CT field-of-view (FOV) | 136 cm (53.3 in)                   |

| SPECT Acquisition                               |                                                                                                   |
|-------------------------------------------------|---------------------------------------------------------------------------------------------------|
| Energy range (photopeak window center)          | 35-588 keV                                                                                        |
| Acquisition modes                               | Static, dynamic, gated, SPECT, gated SPECT, dynamic SPECT, whole-body, whole-body SPECT, SPECT/CT |
| Nuclear Medicine and SPECT Acquisition (Static) |                                                                                                   |
| Time                                            | 50-32000000 ms                                                                                    |

|                                                                                             |                                                                                                               |
|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Counts                                                                                      | 1-2147483647                                                                                                  |
| Zoom                                                                                        | 1.00, 1.23, 1.45, 1.78, 2.00, 2.29, 2.67, 3.20                                                                |
| Matrix                                                                                      | 64x64, 128x128, 256x256, 512x512, 1024x1024                                                                   |
| Body position                                                                               | Supine, prone                                                                                                 |
| Orientation                                                                                 | Head-in, head-out, gurnee-right, gurnee-left, head-left, head-right, sitting, standing, open-right, open-left |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both                                                                              |
| Detector configuration                                                                      | 180°, 90°, 76°, out-facing                                                                                    |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, ME, HE, pinhole                                                                      |

| Nuclear Medicine and SPECT Acquisition (Dynamic)                                            |                                                                                                               |
|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Time                                                                                        | 50-32000000 ms                                                                                                |
| Number of frames                                                                            | 1-2,048 frames                                                                                                |
| Zoom                                                                                        | 1.00, 1.23, 1.45, 1.78, 2.00, 2.29, 2.67, 3.20                                                                |
| Matrix                                                                                      | 64x64, 128x128, 256x256                                                                                       |
| Number of phases                                                                            | 1-32 phases                                                                                                   |
| Body position                                                                               | Supine, prone                                                                                                 |
| Orientation                                                                                 | Head-in, head-out, gurnee-right, gurnee-left, head-left, head-right, sitting, standing, open-right, open-left |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both                                                                              |
| Detector configuration                                                                      | 180°, 90°, 76°, out-facing                                                                                    |
| Acquire with R-wave gate                                                                    | Selectable                                                                                                    |
| Acquire with statics                                                                        | Selectable                                                                                                    |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, ME, HE, pinhole                                                                      |
| Gated                                                                                       |                                                                                                               |
| Time                                                                                        | 1-32000000 ms                                                                                                 |
| Counts                                                                                      | 1-15000000 cts                                                                                                |
| Zoom                                                                                        | 1.00, 1.23, 1.45, 1.78, 2.00, 2.29, 2.67, 3.20                                                                |
| Matrix                                                                                      | 64x64, 128x128                                                                                                |
| Number of frames                                                                            | 2-32 frames                                                                                                   |
| Body position                                                                               | Supine, prone                                                                                                 |

|                                                                                             |                                                                                                               |
|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Orientation                                                                                 | Head-in, head-out, gurnee-right, gurnee-left, head-left, head-right, sitting, standing, open-right, open-left |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both                                                                              |
| Detector configuration                                                                      | 180°, 90°, 76°, out-facing                                                                                    |
| Heartbeats                                                                                  | 1-100000 heartbeats                                                                                           |
| Heart beat framing                                                                          | Forward, forward/backward by thirds                                                                           |
| Beat window % width                                                                         | 0-200                                                                                                         |
| Beat window center                                                                          | 256-2,000 ms/beat                                                                                             |
| Autocenter primary window                                                                   | Selectable                                                                                                    |
| Autotracking                                                                                | Selectable                                                                                                    |
| Reject PVC beats                                                                            | Selectable                                                                                                    |
| Beats to reject post PVC                                                                    | 0-6                                                                                                           |
| PVC threshold (bpm)                                                                         | 1-99 beats                                                                                                    |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, ME, HE, pinhole                                                                      |

|                                                                                             | Whole-body                                     |
|---------------------------------------------------------------------------------------------|------------------------------------------------|
| Scan speed with autocontour                                                                 | 3-60 cm/min                                    |
| Zoom                                                                                        | 1.00                                           |
| Matrix                                                                                      | 256x512, 256x1024, 512x1024                    |
| Scan length                                                                                 | 1-203 cm                                       |
| Body position                                                                               | Supine, prone                                  |
| Orientation                                                                                 | Head out                                       |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both               |
| Detector configuration                                                                      | 180°                                           |
| Autocontour                                                                                 | Selectable                                     |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, ME, HE                |
| SPECT                                                                                       |                                                |
| Time                                                                                        | 500-32000000 ms                                |
| First view by counts                                                                        | 1-100000 kcts                                  |
| Zoom                                                                                        | 1.00, 1.23, 1.45, 1.78, 2.00, 2.29, 2.67, 3.20 |
| Matrix                                                                                      | 64x64, 128x128, 256x256                        |
| Maximum number of views                                                                     | 360 per head                                   |

|                                                                                             |                                                                                    |
|---------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| Body position                                                                               | Supine, prone                                                                      |
| Orientation                                                                                 | Head-in (only 180°), head-out (all configurations)                                 |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both                                                   |
| Detector configuration                                                                      | 180°, 90°, 76°, IQ•SPECT                                                           |
| Orbit                                                                                       | Circular (180°, 90°), NCO (180°, 90°, 76°), NCO-prescan (90°, 76°), cardio-centric |
| Start angle                                                                                 | -179°-180°                                                                         |
| Mode                                                                                        | Step and shoot, continuous, acquire during step                                    |
| Degrees of rotation                                                                         | 90° (only 90°), 104° (only 76°), 180° (90° and 180°), 360° (90° and 180°)          |
| Rotation direction                                                                          | Clockwise, counterclockwise                                                        |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, LEFB, ME, HE, SMARTZOOM                                   |
|                                                                                             | <b>Dynamic SPECT</b>                                                               |
| Time/cycle                                                                                  | 10-900 sec.                                                                        |
| Cycles/repeat                                                                               | 1-10 cycles/repeat                                                                 |
| Repeats/phase                                                                               | 1-80 repeats/phase                                                                 |
| Number of phases                                                                            | 1-16 phases                                                                        |
| Zoom                                                                                        | 1.00, 1.23, 1.45, 1.78, 2.00, 2.29, 2.67, 3.20                                     |
| Matrix                                                                                      | 64x64, 128x128                                                                     |
| Start angle                                                                                 | -179°-180°                                                                         |
| Body position                                                                               | Supine, prone                                                                      |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both                                                   |
| Orientation                                                                                 | Head-in (only 180°), head-out (all configurations)                                 |
| Detector configuration                                                                      | 180°, 90°, 76°, IQ•SPECT                                                           |
| Orbit                                                                                       | Circular (180°, 90°), NCO (180°, 90°, 76°), NCO-prescan (90°, 76°), cardio-centric |
| Mode                                                                                        | Step and shoot, continuous                                                         |
| Rotation direction                                                                          | Clockwise, counterclockwise                                                        |
| Degrees of rotation                                                                         | 90° (only 90°), 104° (only 76°), 180° (90° and 180°), 360° (90° and 180°)          |
| Pause before phase                                                                          | Selectable                                                                         |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, LEFB, ME, HE, SMARTZOOM                                   |

| Whole- body SPECT                                                                           |                                                 |
|---------------------------------------------------------------------------------------------|-------------------------------------------------|
| Time                                                                                        | 500-32000000 ms                                 |
| First view by counts                                                                        | 1-100000 kcts                                   |
| Zoom                                                                                        | 1.00                                            |
| Number of bed positions                                                                     | 2-5 bed positions                               |
| Matrix                                                                                      | 64x64, 128x128, 256x256                         |
| Orientation                                                                                 | Head out                                        |
| Body position                                                                               | Supine, prone                                   |
| Detectors                                                                                   | Detector 3/8, detector 5/8, both                |
| Detector configuration                                                                      | 180°                                            |
| Orbit                                                                                       | Circular (180°, 90°), NCO (180°, 90°, 76°)      |
| Mode                                                                                        | Step and shoot, continuous, acquire during step |
| Rotational direction                                                                        | Clockwise, counterclockwise                     |
| Degrees of rotation                                                                         | 180°, 360°                                      |
| Start angle                                                                                 | -179°-180°                                      |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, ME, HE                 |

| Gated SPECT            |                                                                                    |
|------------------------|------------------------------------------------------------------------------------|
| Time                   | 500-32000000 ms                                                                    |
| Accepted beats/view    | 1-99 beats/view                                                                    |
| Zoom                   | 1.00, 1.23, 1.45, 1.78, 2.00, 2.29, 2.67, 3.20                                     |
| Matrix                 | 64x64, 128x128                                                                     |
| Number of frames       | 2-32 frames                                                                        |
| Body position          | Supine, prone                                                                      |
| Orientation            | Head out                                                                           |
| Detectors              | Detector 3/8, detector 5/8, both                                                   |
| Detector configuration | 180°, 90°, 76°, IQ•SPECT                                                           |
| Orbit                  | Circular (180°, 90°), NCO (180°, 90°, 76°), NCO-prescan (90°, 76°), cardio-centric |
| Start angle            | -179°-180°                                                                         |
| Mode                   | Step and shoot                                                                     |
| Degrees of rotation    | 90° (only 90°), 104° (only 76°), 180° (90° and 180°), 360° (90° and 180°)          |

|                                                                                             |                                            |
|---------------------------------------------------------------------------------------------|--------------------------------------------|
| Rotation direction                                                                          | Clockwise, counterclockwise                |
| Heart beat framing                                                                          | Forward, forward/backward by thirds        |
| Beat window % width                                                                         | 0-200                                      |
| Beat window center                                                                          | 256-2,000 ms/beat                          |
| Autocenter primary window                                                                   | Selectable                                 |
| Autotracking                                                                                | Selectable                                 |
| Reject PVC beats                                                                            | Selectable                                 |
| Beats to reject post PVC                                                                    | 0-6                                        |
| PVC threshold (bpm)                                                                         | 1-99 beats                                 |
| Allowable collimators (all collimators may not be supported by all detector configurations) | LEHR, LPHR, LEAP, LEUHR, ME, HE, SMARTZOOM |

| SPECT Motions                                        |                                |
|------------------------------------------------------|--------------------------------|
| Average Autocontour distance                         | 1.1 cm (0.45 in)               |
| Max. radial and lateral                              | 72 cm/min (28.3 in/min)        |
| Max. lateral position left/right                     | 37.5 cm (14.7 in)/10 cm (4 in) |
| Max. clockwise/counter-clockwise rotation detector 1 | 405°/-135°                     |
| Ring rotation range                                  | 540°                           |
| Rotational uniformity                                | Yes                            |
| Rotational accuracy                                  | 0.1°                           |
| Rotational speed                                     | 0.03-3.0 RPM                   |
| Center of rotation                                   | ≤0.25 pixel (64x64 matrix)     |
| Max. caudal tilt                                     | +16°/-16°                      |

| Patient Bed(table)     |                                |
|------------------------|--------------------------------|
| Width                  | 81.9 cm (32.2 in)              |
| Length                 | 248.0 cm (8 ft 1.6 in)         |
| Weight without ICC/ACC | 950 kg (2,096 lb)              |
| Height                 | 112.0 cm (3 ft 8 in)           |
| Vertical motion range  | 48.0-112.0 cm (19-44 in)       |
| Vertical speed         | 72 cm/min (28 in/min), maximum |
| Pallet material        | Carbon fiber                   |
| Pallet thickness       | 15 mm (.6 in)                  |
| Pallet width           | 40.0 cm (15.8 in)              |

|                                     |                                               |
|-------------------------------------|-----------------------------------------------|
| Attenuation at 140 keV              | <10%                                          |
| Max. patient weight                 | 227 kg (500 lb)                               |
| Max. deflection of patient pallet   | <2.0 mm (<0.08 in) for 92 kg (200 lb) patient |
| Max. scan length in whole-body mode | 203 cm (6 ft 6.7 in)                          |
| Horizontal motion accuracy          | 0.7 mm (0.03 in)                              |
| Min./max. horizontal speed          | 3-600 cm/min (1.2-236 in/min)                 |
|                                     |                                               |
|                                     |                                               |
|                                     |                                               |
| <b>Optional Pallets</b>             |                                               |
| <b>Pediatric</b>                    |                                               |
| Material                            | Carbon fiber composite                        |
| Thickness                           | 0.6 cm (0.25 in)                              |
| Width                               | 25.4 cm (10 in)                               |
| Length                              | 145 cm (57 in)                                |
| Weight                              | 7.3 kg (16 lb)                                |
| Attenuation at 140 keV              | <10%                                          |
| Max. patient weight                 | 27 kg (60 lb)                                 |
| <b>Scintimammography</b>            |                                               |
| Material                            | Carbon fiber composite                        |
| Thickness                           | 1.6 cm (0.63 in)                              |
| Width                               | 35.6 cm (14 in)                               |
| Length                              | 190.5 cm (75 in)                              |
| Weight                              | 7.7 kg (17 lb)                                |
| Attenuation at 140 keV              | <10%                                          |
| Max. patient weight                 | 135 kg (300 lb)                               |
| <b>Optional Pallets</b>             |                                               |
| <b>Radiotherapy planning</b>        |                                               |
| Material                            | Carbon fiber composite                        |
| Thickness                           | 1.5 cm (0.6 in)                               |
| Width                               | 53 cm (20.9 in)                               |
| Length                              | 203.5 cm (80.1 in)                            |
| Weight                              | 9 kg (20 lb)                                  |
| Attenuation at 140 keV              | <10%                                          |

|                               |                 |
|-------------------------------|-----------------|
| Max. patient weight           | 227 kg (500 lb) |
| <b>Pallet holder</b>          |                 |
| Width                         | 58 cm           |
| Length                        | 25 cm           |
| Weight                        | 2 kg            |
| Maximum load bearing capacity | 18.1 kg         |

|                                  |                                           |
|----------------------------------|-------------------------------------------|
| <b>Rear Pallet Support</b>       |                                           |
| Width                            | 26.3 cm (10.3 in)                         |
| Length                           | 104.3 cm (3 ft 5.1 in)                    |
| Weight                           | 188.3 kg (415.2 lb)                       |
| <b>ECG Trigger</b>               |                                           |
| Integration                      | Internal (inside patient bed) or external |
| Framing modes                    | Forward or forward/backward by thirds     |
| Buffered beat window             | Yes                                       |
| Bad beat rejection               | Yes                                       |
| Criteria for framing images      | Frames/R-R interval                       |
| Beat acceptance window           | Automatic or manual selection             |
| <b>Collimator Exchanger Cart</b> |                                           |
| Height                           | 101.4 cm (3 ft 3.9 in)                    |
| Width                            | 82.8 cm (2 ft 8.6 in)                     |
| Depth                            | 120.4 cm (3 ft 11.4 in)                   |
| Weight                           | 181.4 kg (400 lb)                         |
| <b>FOV Dimensions</b>            |                                           |
| FOV                              | 53.3x38.7 cm (21x15.25 in)                |
| Diagonal FOV                     | 65.9 cm (25.9 in)                         |
| <b>Crystal</b>                   |                                           |
| Size                             | 59.1x44.5 cm (23.25x17.5 in)              |
| Diagonal                         | 73.9 cm (29.1 in)                         |
| Thickness                        | 9.5 mm (3/8 in) or 15.9 mm (5/8 in)       |
| <b>Photomultiplier Tubes</b>     |                                           |
| Total number                     | 59                                        |
| Diameter                         | 53-7.6 cm (3 in) and 6-5.1 cm (2.4-2 in)  |

|                                |                                           |
|--------------------------------|-------------------------------------------|
| Type                           | Bialkali high-efficiency box-type dynodes |
| Array                          | Hexagonal                                 |
| Sampling rate                  | 30.0 MHz                                  |
| <b>Detector Shielding</b>      |                                           |
| Back                           | 9.5 mm (0.375 in)                         |
| Sides                          | 12.7 mm (0.5 in)                          |
| Min./max. in patient direction | 27.9/36.4 mm (1.1/1.435 in)               |
| Brain reach <sub>4</sub>       | 7.6 cm (3 in)                             |

| Detector                                                          | 3/8"      | 5/8"      |
|-------------------------------------------------------------------|-----------|-----------|
| Intrinsic spatial resolution                                      |           |           |
| Full width at half maximum (FWHM) in central field of view (CFOV) | ≤3.8 mm   | ≤4.5 mm   |
| FWHM in useful field of view (UFOV)                               | ≤3.9 mm   | ≤4.6 mm   |
| Full width at tenth maximum (FWTM) in CFOV                        | ≤7.5 mm   | ≤8.7 mm   |
| FWTM in UFOV                                                      | ≤7.7 mm   | ≤8.9 mm   |
| Intrinsic spatial linearity                                       |           |           |
| Differential in CFOV                                              | ≤0.2 mm   | ≤0.2 mm   |
| Differential in UFOV                                              | ≤0.2 mm   | ≤0.2 mm   |
| Absolute in CFOV                                                  | ≤0.4 mm   | ≤0.5 mm   |
| Absolute in UFOV                                                  | ≤0.7 mm   | ≤1.0 mm   |
| Intrinsic energy resolution                                       |           |           |
| FWHM in CFOV                                                      | ≤9.9%     | ≤9.9%     |
| Intrinsic flood field uniformity (uncorrected)                    |           |           |
| Differential in CFOV                                              | ≤2.5%     | ≤2.5%     |
| Differential in UFOV                                              | ≤2.7%     | ≤2.7%     |
| Integral in CFOV                                                  | ≤2.9%     | ≤2.9%     |
| Integral in UFOV                                                  | ≤3.7%     | ≤3.7%     |
| Multiple window spatial registration                              | ≤0.6 mm   | ≤1.0 mm   |
| Intrinsic count rate performance in air                           |           |           |
| Maximum count rate                                                | ≥460 kcps | ≥460 kcps |
| Maximum count rate (@15% window)                                  | ≥310 kcps | ≥310 kcps |
| Intrinsic spatial resolution at 75 kcps                           |           |           |
| FWHM in UFOV                                                      | ≤4.1 mm   | ≤4.6 mm   |
| FWTM in UFOV                                                      | ≤7.8 mm   | ≤8.9 mm   |

|                                                           |              |              |
|-----------------------------------------------------------|--------------|--------------|
| Intrinsic flood field uniformity at 75 kcps (uncorrected) |              |              |
| Differential in CFOV                                      | $\leq 2.5\%$ | $\leq 2.5\%$ |
| Differential in UFOV                                      | $\leq 2.7\%$ | $\leq 2.7\%$ |
| Integral in CFOV                                          | $\leq 2.9\%$ | $\leq 2.9\%$ |
| Integral in UFOV                                          | $\leq 3.7\%$ | $\leq 3.7\%$ |

|                                                                                        |                                                                                                 |
|----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| Low-Energy High-Resolution Collimator                                                  |                                                                                                 |
| Type of collimator (ex. Parallel, Focusing, Converging, Diverging, Pinhole)            | Parallel Hex Hole,<br>Hole dia AF: 1.11mm<br>Septal thickness: .16 mm<br>Hole length: 24.05 mm  |
| Over all collimator size Dimensions, mm                                                | Length: 635 mm<br>Width: 531.8 mm<br>Height: 27.2 mm                                            |
| Distance to focus, mm (If collimator is focusing, focal length)                        | N/A non focusing                                                                                |
| Output field of view, mm                                                               | 533 mm x 387 mm                                                                                 |
| Visualization of distribution of radionuclide, keV.<br>Resolution of the energy level. | Using Tc-99 (140keV)<br>Geometric Resolution at 10cm: 6.4mm<br>System Resolution at 10cm: 7.5mm |
| Mass, kg                                                                               | 22.1 kg                                                                                         |
| Low Energy Ultra High Resolution                                                       |                                                                                                 |
| Type of collimator (ex. Parallel, Focusing, Converging, Diverging, Pinhole)            | Parallel Hex Hole,<br>Hole dia AF: 1.16 mm                                                      |

|                                                                                       |                                                                                                  |
|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
|                                                                                       | Septal thickness: .13 mm<br>Hole length: 35.8 mm                                                 |
| Over all collimator size Dimensions, mm                                               | Length: 635 mm<br>Width: 531.8 mm<br>Height: 38.3 mm                                             |
| Distance to focus, mm (If collimator is focusing, focal length)                       | N/A non focusing                                                                                 |
| Output field of view, mm                                                              | 533 mm x 387 mm                                                                                  |
| Visualization of distribution of radionuclide, keV<br>Resolution of the energy level. | Using Tc-99 (140keV)<br>Geometric Resolution at 10cm: 4.6mm<br>System Resolution at 10cm: 6.0 mm |
| Mass, kg                                                                              | 28 Kg                                                                                            |
| <b>Fan beam collimator</b>                                                            |                                                                                                  |
| Type of collimator (ex. Parallel, Focusing, Converging, Diverging, Pinhole)           | Parallel Hex Hole,<br>Hole dia AF: 1.16 mm<br>Septal thickness: .13 mm<br>Hole length: 35.8 mm   |
| Over all collimator size Dimensions, mm                                               | Length: 635 mm<br>Width: 531.8 mm<br>Height: 38.3 mm                                             |

|                                                                                       |                                                                                                              |
|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| Distance to focus, mm (If collimator is focusing, focal length)                       | N/A non focusing                                                                                             |
| Output field of view, mm                                                              | 533 mm x 387 mm                                                                                              |
| Visualization of distribution of radionuclide, keV<br>Resolution of the energy level. | Using Tc-99 (140keV)<br><br>Geometric Resolution at 10cm: 4.6mm<br><br>System Resolution at 10cm: 6.0 mm     |
| Mass, kg                                                                              | 28 Kg                                                                                                        |
| <b>Medium Energy Collimator</b>                                                       |                                                                                                              |
| Type of collimator (ex. Parallel, Focusing, Converging, Diverging, Pinhole)           | Parallel Hex Hole,<br><br>Hole dia AF: 2.94 mm<br><br>Septal thickness: 1.14 mm<br><br>Hole length: 40.64 mm |
| Over all collimator size Dimensions, mm                                               | Length: 635 mm<br><br>Width: 531.8 mm<br><br>Height: 44 mm                                                   |
| Distance to focus, mm (If collimator is focusing, focal length)                       | N/A non focusing                                                                                             |
| Output field of view, mm                                                              | 533 mm x 387 mm                                                                                              |
| Visualization of distribution of radionuclide, keV                                    | Using Ga-67 (300 keV)                                                                                        |

|                                                                                       |                                                                                                     |
|---------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| Resolution of the energy level.                                                       | Geometric Resolution at 10cm: 10.8mm<br>System Resolution at 10cm: 12.5 mm                          |
| Mass, kg                                                                              | 63.5 Kg                                                                                             |
| <b>High energy collimator</b>                                                         |                                                                                                     |
| Type of collimator (ex. Parallel, Focusing, Converging, Diverging, Pinhole)           | Parallel Hex Hole,<br>Hole dia AF: 4 mm<br>Septal thickness: 2 mm<br>Hole length: 59.7 mm           |
| Over all collimator size Dimensions, mm                                               | Length: 635.4 mm<br>Width: 531.8 mm<br>Height: 63.3 mm                                              |
| Distance to focus, mm (If collimator is focusing, focal length)                       | N/A non focusing                                                                                    |
| Output field of view, mm                                                              | 533 mm x 387 mm                                                                                     |
| Visualization of distribution of radionuclide, keV<br>Resolution of the energy level. | Using I-131 (364 keV)<br>Geometric Resolution at 10cm: 13.2mm<br>System Resolution at 10cm: 13.4 mm |
| Mass, kg                                                                              | 124.7 Kg                                                                                            |
| <b>Ultra-high energy collimator</b>                                                   |                                                                                                     |
| Type of collimator (ex. Parallel, Focusing,                                           | Parallel Hex Hole,                                                                                  |

|                                                                                       |                                                                                                     |
|---------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| Converging, Diverging, Pinhole)                                                       | Hole dia AF: 4 mm<br>Septal thickness: 2 mm<br>Hole length: 59.7 mm                                 |
| Over all collimator size Dimensions, mm                                               | Length: 635.4 mm<br>Width: 531.8 mm<br>Height: 63.3 mm                                              |
| Distance to focus, mm (If collimator is focusing, focal length)                       | N/A non focusing                                                                                    |
| Output field of view, mm                                                              | 533 mm x 387 mm                                                                                     |
| Visualization of distribution of radionuclide, keV<br>Resolution of the energy level. | Using I-131 (364 keV)<br>Geometric Resolution at 10cm: 13.2mm<br>System Resolution at 10cm: 13.4 mm |
| Mass, kg                                                                              | 124.7 Kg                                                                                            |
| <b>Pinhole-collimator</b>                                                             |                                                                                                     |
| Type of collimator (ex. Parallel, Focusing, Converging, Diverging, Pinhole)           | Pin Hole, Round hole. Available hole sizes: 4 mm, 6 mm, 8mm                                         |
| Over all collimator size Dimensions, mm                                               | Length: 635 mm<br>Width: 531.8 mm<br>Height: 213.9 mm                                               |
| Distance to focus, mm (If collimator is focusing, focal length)                       | Cone Length: 219.3 mm                                                                               |

|                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Output field of view, mm                                                              | Diameter at base of cone: 220 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Visualization of distribution of radionuclide, keV<br>Resolution of the energy level. | <p>Using Tc-99 (140 keV)</p> <p>Geometric resolution at 10 cm with 4mm aperture: 6.2 mm</p> <p>Geometric resolution at 10 cm with 6mm aperture: 9.3 mm</p> <p>Geometric resolution at 10 cm with 8mm aperture: 12.3 mm</p> <p>System resolution at 10 cm with 4mm aperture: 6.6 mm</p> <p>System resolution at 10 cm with 6mm aperture: 9.5 mm</p> <p>System resolution at 10 cm with 8mm aperture: 12.5 mm</p> <p>Using I-123 (159 keV)</p> <p>Geometric resolution at 10 cm with 4mm aperture: 6.3 mm</p> <p>Geometric resolution at 10 cm with 6mm aperture: 9.3 mm</p> <p>Geometric resolution at 10 cm with 8mm aperture: 12.4 mm</p> <p>System resolution at 10 cm with 4mm aperture: 6.6 mm</p> <p>System resolution at 10 cm with 6mm aperture: 9.5 mm</p> <p>System resolution at 10 cm with 8mm aperture: 12.5 mm</p> <p>Using I-131 (364 keV)</p> <p>Geometric resolution at 10 cm with 4mm aperture: 7.5 mm</p> <p>Geometric resolution at 10 cm with 6mm aperture: 10,6 mm</p> <p>Geometric resolution at 10 cm with 8mm aperture: 13.6</p> |

|          |                                                                                                                                                                                          |
|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|          | mm<br><br>System resolution at 10 cm with 4mm aperture: 7.6 mm<br><br>System resolution at 10 cm with 6mm aperture: 10.7 mm<br><br>System resolution at 10 cm with 8mm aperture: 13.7 mm |
| Mass, kg | 80.3 Kg                                                                                                                                                                                  |
|          |                                                                                                                                                                                          |

| Detector with Collimator                                                | 3/8"        | 5/8"        |
|-------------------------------------------------------------------------|-------------|-------------|
| System spatial resolution without scatter (LEHR at 10 cm)               |             |             |
| FWHM in CFOV                                                            | ≤7.5 mm     | ≤7.8 mm     |
| FWTM in CFOV                                                            | ≤13.6 mm    | ≤14.9 mm    |
| System spatial resolution with scatter (LEHR at 10 cm)                  |             |             |
| FWHM in CFOV                                                            | ≤8.3 mm     | ≤8.9 mm     |
| FWTM in CFOV                                                            | ≤18.6 mm    | ≤19.5 mm    |
| System planar sensitivity (LEHR at 10 cm)                               |             |             |
| Absolute 99mTc                                                          | 202 cpm/μCi | 225 cpm/μCi |
| System planar sensitivity (MELP at 10 cm)                               |             |             |
| Absolute 111In                                                          | 430 cpm/μCi | 565 cpm/μCi |
| Detector with Collimator Tomographic                                    | 3/8"        | 5/8"        |
| Reconstructed spatial resolution without scatter at 15 cm radius (LEHR) |             |             |
| Filtered back projection                                                |             |             |
| Central transaxial                                                      | ≤10.2 mm    | –           |
| Central axial                                                           | ≤10.8 mm    | –           |
| Peripheral radial                                                       | ≤9.8 mm     | –           |
| Peripheral tangential                                                   | ≤8.4 mm     | –           |
| Peripheral axial                                                        | ≤9.0 mm     | –           |
| Reconstructed spatial resolution without scatter at 15 cm radius (LEHR) |             |             |
| Flash 3D iterative reconstruction                                       |             |             |
| Central transaxial                                                      | ≤4.4 mm     | –           |
| Central axial                                                           | ≤4.4 mm     | –           |
| Peripheral radial                                                       | ≤4.0 mm     | –           |

|                                                      |                                         |          |
|------------------------------------------------------|-----------------------------------------|----------|
| Peripheral tangential                                | ≤3.9 mm                                 | –        |
| Peripheral axial                                     | ≤4.2 mm                                 | –        |
| Reconstructed spatial resolution with scatter (LEHR) | Filtered back projection                |          |
| Center                                               | ≤10.7 mm                                | ≤11.5 mm |
| Radial                                               | ≤10.9 mm                                | ≤12.0 mm |
| Tangential                                           | ≤7.9 mm                                 | ≤8.8 mm  |
| Reconstructed spatial resolution with scatter (LEHR) | Flash 3D iterative reconstruction       |          |
| Center                                               | ≤5.8 mm                                 | –        |
| Radial                                               | ≤5.0 mm                                 | –        |
| Tangential                                           | ≤4.1 mm                                 | –        |
| Average volume sensitivity per axial centimeter      |                                         |          |
| LEHR, 99mTc                                          | 12,000 (cts/sec)/(MBq/cm <sup>2</sup> ) | –        |
| Detector sensitivity variation LEHR, 99mTc           | ≤5.0%                                   | –        |

| Detector with Collimator Whole-body Scanning                                                 | 3/8"     | 5/8" |
|----------------------------------------------------------------------------------------------|----------|------|
| Whole-body system spatial resolution without scatter at 10 cm/min scan speed (LEHR at 10 cm) |          |      |
| FWHM perpendicular                                                                           | ≤7.5 mm  | –    |
| FWHM parallel                                                                                | ≤7.9 mm  | –    |
| FWTM perpendicular                                                                           | ≤14.0 mm | –    |
| FWTM parallel                                                                                | ≤14.2 mm | –    |

| Collimators                    | LEHR                       | LPHR                            | LEAP                   | LEUHR                            | LEFB                | MELP                          | HE               | SMART-ZOOM                           |
|--------------------------------|----------------------------|---------------------------------|------------------------|----------------------------------|---------------------|-------------------------------|------------------|--------------------------------------|
|                                | Low Energy High Resolution | Low Penetration High Resolution | Low Energy All Purpose | Low Energy Ultra High Resolution | Low Energy Fan Beam | Medium Energy Low Penetration | High Energy      | IQ•SPECT                             |
| Isotope                        | <sup>99m</sup> Tc          | <sup>123</sup> I                | <sup>99m</sup> Tc      | <sup>99m</sup> Tc                | <sup>99m</sup> Tc   | <sup>67</sup> Ga              | <sup>131</sup> I | <sup>99m</sup> Tc                    |
| Hole shape                     | Hex                        | Hex                             | Hex                    | Hex                              | Hex                 | Hex                           | Hex              | Hex                                  |
| Number of holes (x1000)        | 148                        | 86                              | 90                     | 146                              | 64                  | 14                            | 8                | 48                                   |
| Hole length                    | 24.05 mm                   | 35.0 mm                         | 24.05 mm               | 35.8 mm                          | 35 mm               | 40.64 mm                      | 59.7 mm          | 40.25 mm                             |
| Septal thickness               | 0.16 mm                    | 0.2 mm                          | 0.2 mm                 | 0.13 mm                          | 0.16 mm             | 1.14 mm                       | 2 mm             | 0.2-0.4                              |
| Hole diameter across the flats | 1.11 mm                    | 1.5 mm                          | 1.45 mm                | 1.16 mm                          | 1.53 mm             | 2.94 mm                       | 4 mm             | 1.9 mm                               |
| Sensitivity at 10 cm           | 202 cpm/μCi                | 330 cpm/μCi                     | 330 cpm/μCi            | 100 cpm/μCi                      | 280 cpm/μCi         | 275 cpm/μCi                   | 135 cpm/μCi      | 285 cpm/μCi<br>810 cpm/ μCi at 28 cm |
| Geometric resolution at 10 cm  | 6.4 mm                     | 6.4 mm                          | 8.3 mm                 | 4.6 mm                           | 6.3 mm              | 10.8 mm                       | 13.2 mm          | 6.95 mm                              |

| Collimators                | LEHR              | LPHR            | LEAP              | LEUHR           | LEFB              | MELP               | HE                | SMART-ZOOM       |
|----------------------------|-------------------|-----------------|-------------------|-----------------|-------------------|--------------------|-------------------|------------------|
| System resolution at 10 cm | 7.5 mm            | 8.0 mm          | 9.4 mm            | 6.0 mm          | 7.3 mm            | 12.5 mm            | 13.4 mm           | 7.4 mm8          |
| Calculated penetration     | 1.5%              | 1.2%            | 1.9%              | 0.8%            | 1.0%              | 1.2%               | 3.5%              | N/A              |
| Weight                     | 22.1 kg (48.7 lb) | 33.1 kg (73 lb) | 22.6 kg (49.8 lb) | 28 kg (61.8 lb) | 28.4 kg (62.5 lb) | 63.5 kg (140.1 lb) | 124.7 kg (275 lb) | 47.2 kg (104 lb) |

| Pinhole Collimator                      | Isotope              |                      |                      |
|-----------------------------------------|----------------------|----------------------|----------------------|
|                                         | 99mTc                | 123I                 | 131I                 |
| Hole shape                              | Round                | Round                | Round                |
| Number of holes                         | 1                    | 1                    | 1                    |
| Cone aperture                           | 4 mm<br>6 mm<br>8 mm | 4 mm<br>6 mm<br>8 mm | 4 mm<br>6 mm<br>8 mm |
| Cone length                             | 219.3 mm             | 219.3 mm             | 219.3 mm             |
| Diameter at base of cone (approximate)  | 220 mm               | 220 mm               | 220 mm               |
| Sensitivity at 10 cm with 4 mm          | 123 cpm/ $\mu$ Ci    | 111 cpm/ $\mu$ Ci    | 67 cpm/ $\mu$ Ci     |
| Sensitivity at 10 cm with 6 mm          | 271 cpm/ $\mu$ Ci    | 243 cpm/ $\mu$ Ci    | 133 cpm/ $\mu$ Ci    |
| Sensitivity at 10 cm with 8 mm          | 478 cpm/ $\mu$ Ci    | 426 cpm/ $\mu$ Ci    | 221 cpm/ $\mu$ Ci    |
| Geometric resolution at 10 cm with 4 mm | 6.2 mm               | 6.3 mm               | 7.5 mm               |
| Geometric resolution at 10 cm with 6 mm | 9.3 mm               | 9.3 mm               | 10.6 mm              |
| Geometric resolution at 10 cm with 8 mm | 12.3 mm              | 12.4 mm              | 13.6 mm              |
| System resolution at 10 cm with 4 mm    | 6.6 mm               | 6.6 mm               | 7.6 mm               |
| System resolution at 10 cm with 6 mm    | 9.5 mm               | 9.5 mm               | 10.7 mm              |
| System resolution at 10 cm with 8 mm    | 12.5 mm              | 12.5 mm              | 13.7 mm              |
| Weight                                  | 80.3 kg (177 lb)     | 80.3 kg (177 lb)     | 80.3 kg (177 lb)     |

| Gantry Dimensions                          |                                  |
|--------------------------------------------|----------------------------------|
| Aperture                                   | 70 cm                            |
| Scan field                                 | 50 cm                            |
| Rotation time                              | 0.6 s<br>0.8 s<br>1.0 s<br>1.5 s |
| Temporal resolution (min.)                 | 150 ms                           |
| Data Acquisition System                    |                                  |
| Max. number of slices/rotation             | 6                                |
| Number of physical detector rows           | 16                               |
| Number of physical detector channels/slice | 736                              |

|                                        |                                                |
|----------------------------------------|------------------------------------------------|
| Number of detector elements            | 11,776                                         |
| Total channels per slice               | 1,472                                          |
| Number of projections                  | Up to 1,875<br>(1 s/360°)                      |
| Sequence acquisition modes             | 6x1 mm 6x2 mm 6x3 mm 2x5 mm                    |
| Spiral acquisition modes               | 6x0.5 mm 1x1 mm 6x1 mm 6x2 mm 6x3 mm<br>2x5 mm |
| <b>Tube Assembly</b>                   |                                                |
| Tube                                   | DURA 422<br>MVHigh performance CT X-ray tube   |
| Tube current                           | 20-345 mAs                                     |
| Tube voltage                           | 80, 110, 130 kV                                |
| Tube anode heat storage capacity       | 5.0 MHU                                        |
| Focal spot size according to IEC 60336 | 0.8x0.5 mm/7°<br>0.8x0.7 mm/7°                 |
| <b>Laser</b>                           |                                                |
| Wave length                            | 650 nm                                         |
| Maximum power                          | 0.39 mW                                        |
| Pulse pattern                          | continuous wave                                |

|                    |                                   |
|--------------------|-----------------------------------|
| <b>CARE Filter</b> |                                   |
| CARE filter tube   | Equivalent to 5.5 mm Al at 140 kV |

|                                  |           |
|----------------------------------|-----------|
| CARE filter beam limiting device | 0.5 mm Al |
| <b>Generator</b>                 |           |
| Max. power                       | 50 kW     |

|                                         |                                                                            |
|-----------------------------------------|----------------------------------------------------------------------------|
| <b>Topogram</b>                         |                                                                            |
| Length (max.)                           | 184 cm (6 ft)                                                              |
| Scan times                              | 2.1-19.3 s                                                                 |
| Views                                   | a.p., p.a., lateral                                                        |
| <b>Sequence Acquisition</b>             |                                                                            |
| Reconstructed slice widths              | 1.0, 2.0, 3.0,<br>4.0, 5.0, 6.0,<br>9.0, 10.0, 12.0,<br>18.0 mm            |
| Scan times full scan (360°)             | 0.6, 0.8, 1.0,<br>1.5 s (±5%)                                              |
| Partial scan times (240°)               | 0.4, 0.53 s (±5%)                                                          |
| Number of uninterrupted scans per range | 99                                                                         |
| Number of ranges in autorange           | 8                                                                          |
| Standard scan cycle time (±10%)         | 2.1 s at 0.6 s scan time,<br>2.4 s at 0.8 s scan time                      |
| <b>Multislice Spiral Acquisition</b>    |                                                                            |
| Reconstructed slice widths              | 0.63, 0.75, 1.0,<br>1.25, 2.0, 2.5,<br>3.0, 4.0, 5.0, 6.0,<br>8.0, 10.0 mm |
| Scan times full scan (360°)             | 0.6, 0.8, 1.0, 1.5 s                                                       |
| Reconstruction increment                | 0.1-10 mm                                                                  |
| Pitch factor                            | 0.4-1.8                                                                    |

|                                                   |                                                    |
|---------------------------------------------------|----------------------------------------------------|
| Volume pitch                                      | 3.0-10.8                                           |
| Spiral scan time max.                             | 100 s                                              |
| CT scan range                                     | 0-200 cm                                           |
| Continuous scan length and SPECT/CT co-scan range | 168 cm (5 ft 6 in)                                 |
| <b>Dynamic Multiscan</b>                          |                                                    |
| Dynamic scan cycle time ( $\pm 10\%$ )            | 0.9 s at 0.6 s scan time, 1.2 s at 0.8 s scan time |

|                                 |                       |
|---------------------------------|-----------------------|
| <b>Image Reconstruction</b>     |                       |
| Real-time display               | 512x512               |
| Slice thickness                 | 1.0-18.0 mm           |
| Scan field                      | 50 cm                 |
| Recon field                     | 5-50 cm,<br>5-70 cm   |
| Recon time                      | up to 8 images/s      |
| Recon matrix                    | 512x512               |
| HU scale                        | -1,024 to +3,071      |
| Extended HU scale               | -10,240 to<br>+30,710 |
| <b>Phantom CATPHAN (16 cm)</b>  |                       |
| Object size                     | 3 mm                  |
| Contrast difference             | 3 HU                  |
| Dose at surface                 | 19.7 mGy at 100 mAs   |
| Technique                       | 0.6 s, 10 mm, 130 kV  |
| <b>Phantom CATPHAN (20 cm)</b>  |                       |
| Object size                     | 5 mm                  |
| Contrast difference             | 3 HU                  |
| Dose at surface                 | 15.8 mGy at 90 mAs    |
| Technique                       | 0.6 s, 10 mm, 130 kV  |
| <b>High-contrast Resolution</b> |                       |
| 50% MTF ( $\pm 10\%$ )          | 17.5 lp/cm,           |

|                                                                                        |                                         |
|----------------------------------------------------------------------------------------|-----------------------------------------|
|                                                                                        | 0.29 mm                                 |
| 2% MTF ( $\pm 10\%$ )                                                                  | 15.1 lp/cm,<br>0.32 mm                  |
| Technique: Tungsten wire in air                                                        | 160 mAs, 130 kV,<br>0.8 s, 1.0 mm       |
| <b>Homogeneity</b>                                                                     |                                         |
| Cross-field uniformity in a 20 cm water phantom positioned near the center of rotation | Typical $\pm 2$ HU<br>(max. $\pm 4$ HU) |
| <b>Dose, CTD<sub>100</sub> Values</b>                                                  |                                         |
| Phantom $\varnothing$                                                                  |                                         |
|                                                                                        | 110 kV 130 kV                           |
| 16 cm A                                                                                | 13.3 20.1                               |
| B                                                                                      | 13.6 20.3                               |
| 32 cm A                                                                                | 3.9 6.2                                 |
| B                                                                                      | 7.6 11.6                                |

|                                                        |                                   |
|--------------------------------------------------------|-----------------------------------|
| <b>MI Workplace acquisition console</b>                |                                   |
| Processor                                              | Intel® Xeon® Processor E5-1620 v4 |
| Clock frequency of processor, MHz                      | 3.50 GHz                          |
| The speed of information processing                    | 2400 MHz                          |
| Number of cores                                        | 4                                 |
| Hard disk                                              | 1 TB                              |
| Value of memory, GB                                    | 16 GB                             |
| <b>Monitor of the MI Workplace acquisition console</b> |                                   |
| Brightness, nit                                        | 250 cd/m <sup>2</sup>             |
| Contrast range                                         | 1000:1 typ.                       |
| Resolution                                             | 1280 x 1024                       |
| Diagonal, mm                                           | 481 mm                            |

|                                                                      |                                      |
|----------------------------------------------------------------------|--------------------------------------|
| Mass, kg                                                             | 5.9 kg                               |
| Dimensions, mm                                                       | 405 mm (W) x 506 mm (H) x 205 mm (D) |
| <b>Monitor, PPM</b>                                                  |                                      |
| Resolution                                                           | 1280x1024                            |
| Diagonal, mm                                                         | 431.8                                |
| Mass, kg                                                             | 5.0                                  |
| Dimensions, mm                                                       | 380 mm (W) x 353 mm (H) x 58 mm (D)  |
| <b>MultiModality Workplace</b>                                       |                                      |
| Processor                                                            | Intel® Xeon® Processor E3-1275 v3    |
| Clock frequency of processor, MHz                                    | 3.60 GHz                             |
| The speed of information processing                                  | 2133 MHz                             |
| Number of cores                                                      | 4                                    |
| Hard disk                                                            | 1 TB                                 |
| Value of memory, GB                                                  | 8 GB                                 |
| Monitor                                                              | -                                    |
| Resolution                                                           | 1280 x 1024                          |
| Diagonal, mm                                                         | 480 mm                               |
| Soft ware                                                            | -                                    |
| Version                                                              | VB20                                 |
| Production date                                                      | 03/05/2018                           |
| Class of safety accordance to IEC 62304                              | B                                    |
| <b>Symbia.net (PC hardware, monitor, keyboard, mouse, software).</b> |                                      |
| Main block                                                           | -                                    |

|                                          |                                                                             |
|------------------------------------------|-----------------------------------------------------------------------------|
| Mass, kg                                 | 25 kg                                                                       |
| Dimensions, mm                           | 203.2 mm (W) x 444.5 mm (H) x 525.8 mm (D)                                  |
| Processor                                | Intel® Xeon® Processor E5-2637 v4, dual CPU                                 |
| Clock frequency of processor, MHz        | 3.5 GHz                                                                     |
| The speed of information processing      | 2400 MHz                                                                    |
| Number of cores                          | 4 cores,                                                                    |
| Hard disk                                | 480 GB                                                                      |
| Value of memory, GB                      | 32 GB                                                                       |
| Monitor                                  | -                                                                           |
| Resolution                               | 1280 X1024                                                                  |
| Diagonal, mm                             | 480 mm                                                                      |
| Soft ware                                | -                                                                           |
| Version                                  | VB20                                                                        |
| Production date                          | 03/05/2018                                                                  |
| Class of safety accordance to IEC 62304  | B                                                                           |
| <b>Rectangular Fillable Flood Source</b> |                                                                             |
| Dimensions, mm                           | Dimensions 20½" x 28" x 1¼" thick<br>Cavity Dimensions 18½" x 26" x ½" deep |
| Mass, kg                                 | 9.07 kg                                                                     |
| <b>UPS</b>                               |                                                                             |
| UPS (Gantry)                             |                                                                             |
| Nominal power, VA/W                      | 5KVA/3500W                                                                  |

|                     |                                      |
|---------------------|--------------------------------------|
| Input voltage, V    | 208VAC                               |
| Dimensions, mm      | 432, 130, 660 HxWxD                  |
| Mass, kg            | 55                                   |
| UPS (computers)     |                                      |
| Nominal power, VA/W | 2200VA / 1500W                       |
| Input voltage, V    | 230Vac                               |
| Dimensions, mm      | 440 mm (W) x 131 mm (H) x 490 mm (D) |
| Mass, kg            | 24.5 kg                              |
|                     |                                      |

### 8.3. Norms of Radiation Safety

The Symbia T6 combined single photon emission and X-ray computed tomography system with accessories and its operation comply with the radiation safety standards of the Radiation Safety Standards (HPB-99/2009) SanPiN 2.6.1.2523-09.

This product may contain radioactive material or create ionizing radiation.

Ionizing radiation is a hazard for both the operator and all those who are near equipment if radiation precautions are not taken.

To avoid excessive exposure to radiation, clean the room for scanning from the personnel during maintenance and associated scanning. Adhere to the procedures used in your medical institution.

Work and scattered radiation can cause serious or fatal injuries to patients and others if the equipment is used by an unqualified operator. Sufficient precautions should always be observed to avoid exposing the beam of working radiation, as well as indirect radiation, including scattered radiation from the inner parts of the product and any obstacles lying on the path of the beam.

Persons authorized to operate the equipment and directly working with the system or exercising control over its operation must carefully read and strictly observe the established safe dosage standards of exposure, as well as familiarize themselves with the procedures described in the literature,

Failure to observe these precautions can result in death or serious injury to the operator or personnel near the equipment.

#### Precautions associated with ionizing radiation

Radiation safety of personnel and the public from sources of potential exposure is provided by the application of technical measures to reduce the probability of events, due to which the boundary values of the generalized risk can be exceeded.

To protect the operator and the patient being examined, the following precautions against ionizing radiation must be taken.

The doctor performing the examination is responsible for protecting the patient from excessive ionizing radiation.

If necessary, use protective screens to minimize the dose of irradiation of unprotected organs.

Always follow local safety guidelines.

To protect personnel working with a medical device, the following must be ensured:

- receipt of a sanitary and epidemiological certificate for the products that contain ionizing radiation sources;
- development of control levels of radiation factors in the organization and sanitary protection zone, as well as instructions for radiation safety and instructions on the actions of personnel in radiation accidents;
- a list of persons belonging to the personnel has been established;
- conditions for working with sources of ionizing radiation that comply with the rules are created;
- planning and implementation of measures to ensure and improve radiation safety in the organization;
- systematic control of the radiation situation at workplaces, in the premises, on the territory of the organization, in the sanitary protection zone and in the observation zone, as well as the release and discharge of radioactive substances;
- control and accounting of individual doses to personnel;
- informing staff about radiation levels at workplaces and individual doses of radiation;
- training and certification on the issues of ensuring radiation safety of managers and executors of works, specialists of radiation safety services, other persons who constantly or temporarily perform work with sources of ionizing radiation;
- conduct of briefing and testing of personnel's knowledge in the field of radiation safety;
- Conducting preliminary (upon admission to work) and periodic medical examinations of personnel;

- an annual submission of the completed radiation and hygiene passport of the organization in a timely manner;

Personnel working with this medical product should know and strictly comply with radiation safety requirements established by sanitary norms and regulations in force in the territory of the Russian Federation:

- use personal protective equipment in specified cases;
- comply with the established requirements for the prevention of a radiation accident and the rules of conduct in the event of its occurrence;
- timely pass periodic medical examinations and follow the recommendations of the medical commission;
- immediately notify the head (workshop, site, laboratory) and the radiation safety service (the person responsible for radiation safety) of all detected malfunctions in the operation of plants, instruments and apparatus that are sources of ionizing radiation;

When using the product on the territory of the Russian Federation, the annual effective dose of personnel exposure should not exceed 20 mSv per year on average for any consecutive 5 years, but not more than 50 mSv per year.

For women under the age of 45 years working with radiation sources, additional restrictions are introduced: the equivalent dose on the surface of the lower part of the abdominal region should not exceed 1 mSv per month.

For students and students over 16 years of age who are in vocational training using radiation sources, annual doses should not exceed 5 mSv per year on average for any consecutive 5 years, but not more than 12.5 mSv per year.

## 9. Capabilities of the Medical Device with other medical devices:

To ensure the necessary safety, only products or components expressly approved by Siemens, may be used in combination with the system. The person connecting additional equipment to the medical device is considered to be the system configurator and is therefore responsible for ensuring that the current system configuration complies with the relevant standards (for example system standard IEC/EN 60601-ff and/or other applicable standards). In case of doubt, please consult your local contact person.

Example:

Accessories connected to the analog or digital interfaces must be certified according to the respective IEC standards. IEC 60950-1 for Information Technology Equipment and IEC 60601-1 for Medical Equipment. Furthermore all configurations shall comply with the valid version of the system standard IEC 60601-1. Everybody who connects additional equipment to the signal input part configures a medical system. This person is thus responsible that the system complies with the requirements of the valid version of the system standard IEC 60601-1.

If in doubt, consult the Customer Services department or your local representative.


## 10. Sterility:


The medical device is non-sterile.



## 11. Labeling information:

| Symbol                                                                              | Description         |
|-------------------------------------------------------------------------------------|---------------------|
|  | Off and On Symbols  |
|  | Dangerous Voltage   |
|  | Radiation Exposure  |
|  | Do Not Open         |
|  | Alternating Current |
|  | Protective Earth    |

|                                                                                                                                                                                                                                                                                      |                                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|
|                                                                                                                                                                                                     | Emergency Stop                            |
|                                                                                                                                                                                                     | Potential Hand Injur                      |
|                                                                                                                                                                                                     | Type B Applied Part                       |
|                                                                                                                                                                                                     | Defibrillation Proof Type BF Applied Part |
| 1                                                                                                                                                                                                                                                                                    | Detector Designation 3/8                  |
| 2                                                                                                                                                                                                                                                                                    | Detector Designation 5/8                  |
|                                                                                                                                                                                                   | Patient Bed Weight Rating                 |
|                                                                                                                                                                                                   | Whole Body Armrests                       |
| <br><b>CAUTION</b> <br><b>Radioactive Material</b><br>Line Source Gd-153 10 mCi<br>Point Source Co-57 0.05 mCi | Radioactive Material Caution              |


|                                                                                     |                                                                         |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
|    | X-Ray Warning                                                           |
|   | CT System Energized Warning Label                                       |
|  | Laser Radiation                                                         |
|  | This symbol indicates that the cable connection is for an external ECG. |
|  | SNAC Power Switch                                                       |
|  | Camera Power Switch                                                     |

|                                                                                     |                        |
|-------------------------------------------------------------------------------------|------------------------|
|    | Motion Power Switch    |
|    | Do Not Lean Label      |
|   | ECG Label              |
|  | Caution Label Pallet   |
|  | Temperature limitation |
|  | Humidity limitation    |

|                                                                                     |                                         |
|-------------------------------------------------------------------------------------|-----------------------------------------|
|    | Pile legs number limitation             |
|    | Fragile, please handle with care        |
|    | Keep dry                                |
|   | Top. Indicates correct upright position |
|  | Recycling material                      |

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## 12. Transportation and Storage conditions and requirements:

| Parameter           | Range value                   |
|---------------------|-------------------------------|
| Temperature         | -40F (-20°C) to 1200F (+49°C) |
| Temperature Change  | 70F/hour (40C/hour)           |
| Relative humidity   | 10% to 90 % relative humidity |
| Barometric pressure | 700 -1060 hPa                 |

## 13. Service life.

The service life of the system is 10 years.

## 14. Procedure of recycling and disposal.

### 14.1. Requirements for the protection of the environment in the application of medical device

This device does not harm the environment if properly operated, shipped and stored. To maintain compliance with the declared period of harmless use, operation of the equipment must be carried out in strict accordance with the instructions and environmental conditions as described in the operating manual, as well as strict observance of periodic maintenance schedules specified in the "Maintenance procedures" section. The environmentally harmless lifetime for consumables or individual parts may be less than the value established for the device. Periodic replacement of these consumables or parts in order to ensure an environmentally friendly device operation during declared safety period must be performed in accordance with the maintenance rules.

### 14.2. Utilization method/ Disposal

The system must be disposed of in accordance to local regulations for disposal of such products. If you want to stop using the product, please note that municipal public order may contain special rules regarding the disposal of this equipment. To ensure compliance with these legal provisions and to avoid the environmental damage that can be done when you dispose of your system, contact Siemens technical support. Packaging material must be disposed of, taking into account the ecological requirements of

national legislation. From systems and accessories, you need to get rid of at the end of their service life in accordance with applicable regional, national and local regulations:

- Accumulators and batteries
- Capacitors
- Flat Monitors
- Phantoms

For more information, you must contact your local service representative or regional office of Siemens.

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <p>The use of this symbol indicates that this product should not be treated as household waste. By ensuring that this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste-handling of this product. For more detailed information concerning the return and recycling of this product, please consult the supplier from whom you purchased the product</p> |
|  | <p>Concerning Batteries (Directive 2006/66/EC) The directive 2006/66/EC requires separate collection and appropriate disposal of spent batteries. This product contains batteries that are not intended to be replaced by the user. Replacement of those batteries will usually be done during regular maintenance or service by service staff who can also arrange proper disposal.</p>                                                                                                    |

Disposal of the system components, posing a potential danger to human health or the environment, should be carried out with the necessary precautions, and in accordance with legal regulations in force in the territory of the Russian Federation.

For more information on the return and disposal of products and their components or accessories, please contact the local service center or regional office Siemens.

#### *Dispose of packaging materials*

Materials used in packaging equipment, suitable for processing. They should be collected and handled in accordance with the regulations in force in the territory of the Russian Federation.

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## 15. Warranty.

The warranty period for Symbia T6 combined single photon emission and X-ray computed tomography system with accessories is within 12 months from the date of putting the equipment into operation, but not more than 15 months from the date of its shipment. The warranty extends on any malfunctions of medical equipment occurred through the manufactures fault. Repair or replacement of faulty accessories during the warranty period is carried out free of charge.

## 16. Contacts for reclamations.

In case of any questions or problems with handling the system, please use the following contacts:

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC  
Dubinskaya street 96  
115093 Moscow,  
Russia  
Tel.: +7 495 737 10 00  
Fax: +7 495 737 13 20  
[www.healthcare.siemens.ru](http://www.healthcare.siemens.ru)

Contacts for questions and reclamations in the country of Legal manufacturer:

Siemens Medical Solutions USA, Inc.  
Molecular Imaging  
2501 N. Barrington Road,  
Hoffman Estates, Illinois 60192  
USA  
Tel.: +1 847 304 7700  
Email: [www.siemens.com/healthcare](mailto:www.siemens.com/healthcare)

## 17. List of documentation supplemented by this appendix.

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