Reagent kit PraxisTM Extended RAS Panel for qualitative determination of the status of specific mutations genes of RAS (KRAS (exons 2, 3, and 4) and NRAS (exons 2, 3, and 4) by sequencing a new generation in a human genomic DNA sample isolated from colorectal cancer tissue samples fixed in formalin and enclosed in paraffin blocks

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

personally appeared Anaelia Brownridge,
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are
subscribed to the within instrument and acknowledged to me that he/she/they executed the same in
his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the
person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

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Praxis® Extended RAS Panel Reference Guide

FOR IN VITRO DIAGNOSTIC USE

Introduction	3
Getting Started	6
Assay Workflow	15
Qualification of FFPE DNA Samples	16
Set Parameters	21
Hybridization of Oligonucleotide Pool	23
Removal of Unbound Oligonucleotides	26
Extension-Ligation of Bound Oligonucleotides	29
PCR Amplification	30
Verify Library Preparation	34
PCR Clean-Up	36
Library Normalization	39
Library Pooling	42
Library Sequencing	46
Technical Assistance	49

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This product may be used as a component of human molecular diagnostic assays where applicable country laws allow. This product by itself does not provide any diagnostic result. This product includes GoTaq[®] Hot Start Polymerase manufactured by Promega Corporation for distribution by Illumina, Inc. Licensed to Promega Corporation under U.S. Patent Nos. 5,338,671 and 5,587,287 and their corresponding foreign patents.

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
The Praxis™ Extended RAS Panel is a qualitative *in vitro* diagnostic test using targeted high throughput parallel sequencing for the detection of 56 specific mutations in RAS genes [KRAS (exons 2, 3, and 4) and NRAS (exons 2, 3, and 4)] in DNA extracted from formalin-fixed, paraffin-embedded (FFPE) colorectal cancer (CRC) tissue samples. The test is intended to be used on the Illumina MiSeqDx® instrument.

This reference guide provides more detailed instructions, technique tips, and helpful hints to newly trained users to help guide proper execution of the Praxis Extended RAS Panel protocol. This is intended as a supplement and is not meant to replace the package insert.

Pairs of oligonucleotides targeting specific mutations in RAS genes are hybridized to genomic DNA. Each genomic DNA sample is tested with two oligonucleotide pools that independently target each strand of DNA followed by extension and ligation to form DNA templates consisting of the regions of interest flanked by universal primer sequences. Using indexed primers supplied with the kit, DNA templates are then amplified by PCR, pooled into a single tube, and sequenced on the MiSeqDx instrument.

Process Overview

The Illumina Praxis Extended RAS Panel process can be summarized into the following steps:

Create a Run

First, create a run. The MiSeqDx uses the run information to identify each sample and its corresponding indexes.

To create a run, use the Local Run Manager and the Praxis Extended RAS Panel analysis module. Enter the sample name and indexes for the 96-well plate in the module. You can print out the plate layout from the module and use the information for setting up the plate during the assay workflow.

Prepare Libraries

Prepare the libraries using the protocol detailed in this user guide.

Sequence Samples on the MiSeqDx

The Praxis Extended RAS Panel must be sequenced on a MiSeqDx instrument using a paired-end 141 cycle run with dual indexing. For instructions for performing a sequencing run on the MiSeqDx, see the *MiSeqDx Instrument Reference Guide* (document # 1000000021961).

Automated Sequencing and Data Analysis

The Praxis Extended RAS Panel analysis module processes the base calls generated by the MiSeqDx instrument. It is an on-instrument software, which is built in to the instrument's processes. The analysis module produces details indicating whether each sample contains any of the Praxis Extended RAS Panel mutations. For more information about this software, see the *Local Run Manager Software Reference Guide for MiSeqDx* (document # 1000000011880) and the *Local Run Manager Extended RAS Panel Analysis Module Guide* (document # 1000000011875).

Tracking Tools

Illumina provides the following tools for sample tracking and guidance in the lab:

- The **Lab Tracking Form** can be used to record information such as operator name, sample and index information, start and stop times, reagent lot numbers, and barcodes.
- The **Praxis Extended RAS Panel** analysis module is used to create a run. The analysis module provides a feature for recording parameters for the sample plate, such as sample name and dual indexes applicable to the run.

Getting Started

This section describes the Illumina Praxis Extended RAS Panel kit contents, consumables and equipment used, DNA input recommendations, and best practices to apply during the protocol.

Praxis Extended RAS Panel Kit Contents

The Illumina Praxis Extended RAS Panel Kit (Catalog # 20012431) contains the following components. Store the kit components at the specified temperature and in designated pre- and post-amplification areas.

Because the pre- and post-amp reagents are shipped together, it is important to unpack the reagents in the pre-amp lab area, and then move the post-amp reagents to the proper post-amp storage area.

Praxis Extended RAS Panel Box

Box 1A Pre-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Hybridization Buffer	1 tube	4.32 ml	Buffered aqueous solution containing salts and formamide	-25°C to -15°C
Extension-Ligation Mix	1 tube	4.8 ml	Buffered aqueous solution containing proprietary blend of DNA polymerases, DNA ligase, and dNTPs	-25°C to -15°C
Index Primers A-H (A501-A508)	1 tube per primer	192 µl	PCR primers with index sequences and sequencing adapters	-25°C to -15°C
Index Primers 1-12 (A701-A712)	1 tube per primer	128 µl	PCR primers with index sequences and sequencing adapters	-25°C to -15°C
PCR Polymerase	1 tube	56 µl	Proprietary DNA polymerase	-25°C to -15°C
PCR Master Mix	1 tube	2.8 ml	Buffered aqueous solution containing salts and dNTPs	-25°C to -15°C

Box 1B Post-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Library Normalization Diluent	1 tube	4.6 ml	Buffered aqueous solution containing salts, 2-Mercaptoethanol, and formamide	-25°C to -15°C
Library Dilution Buffer	2 tubes	4.5 ml	Buffered aqueous solution	-25°C to -15°C
PhiX Internal Control	1 tube	10 µl	Buffered aqueous solution containing PhiX genomic DNA	-25°C to -15°C

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Praxis Extended RAS Panel, Box 6

Box 6 Pre-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Oligo Pool A–Extended RAS	1 tube	150 µl	Buffered aqueous solution containing oligonucleotides targeting the RAS genes	-25°C to -15°C
Oligo Pool B–Extended RAS	1 tube	150 µl	Buffered aqueous solution containing oligonucleotides targeting the RAS genes	-25°C to -15°C
Quality Control (QC) Primers	4 tubes	75 µl	Buffered aqueous solution containing oligonucleotides for DNA sample qualification	-25°C to -15°C
Quality Control (QC) Template	1 tube	220 µl	Buffered aqueous solution containing genomic DNA for sample qualification	-25°C to -15°C
qPCR Master Mix	2 tubes	1 ml	Buffered aqueous solution containing qPCR reaction components	-25°C to -15°C
Panel Control–Extended RAS	1 tube	8 µl	Buffered aqueous solution containing genomic DNA for assay positive control	-25°C to -15°C

- 10 N NaOH (prepare from tablets or use a standard solution)
- TE Buffer
- RNase/DNase-free water

Post-Amp Reagents

- 10 N NaOH (prepare from tablets or use a standard solution)
- Ethanol, 200 proof for molecular biology
- TE Buffer
- RNase/DNase-free water

Equipment and Materials

Equipment and Materials Required Subsequent to

- 1 **MiSeqDx Instrument**, Catalog # DX-410-1001
- 2 **TruSeq® Index Plate Fixture Kit**, Catalog # FC-130-1005
- 3 **TruSeq Index Plate Fixture & Collar Kit**, Catalog # FC-130-1007

Pre-Amp Equipment and Materials

- 1 **DNA Extraction and Purification Reagents**—See the package insert for more information about DNA Extraction.
- 2 **Heat Block**—One heat block for a 96-well PCR plate is required. The heat block must meet the following performance specifications. Heat blocks with heated lids are required for use.
 - ✦ Temperature Range: Ambient +5°C to 99°C
 - ✦ Temperature Regulation: ±0.1°C at 37°C; ±0.4°C at 60°C
- 3 **Sample Incubator**—One incubator (hybridization oven) is required. The incubator must meet the following performance specifications.
 - ✦ Temperature Range: 10°C to 100°C
 - ✦ Temperature Regulation: ±0.2°C
- 4 **Tabletop Centrifuge**—One tabletop centrifuge is required. (A separate centrifuge is required in the post-amp area.) Any plate centrifuge that attains the designated speeds of the protocol (280 to 2400 × g) is acceptable.
- 5 **Microcentrifuge**—One microcentrifuge for microcentrifuge tubes.
- 6 **Precision Pipettes**—One set of precision pipettes is required. (A separate set is required in the post-amp area.) The use of precision pipettes is required to ensure accurate reagent and sample delivery. If they are calibrated regularly and are accurate within 5% of stated volume, single-channel or multi-channel pipettes can be used.
- 7 **Sealing Rollers or Sealing Wedge**—A sealing roller or sealing wedge are required for sealing microplates.
- 8 **Consumables**—The following consumables are required.
 - ✦ 96-well or 384-well optical PCR plates and seals
 - ✦ 96-well skirted PCR plates, 0.2 ml, polypropylene, or equivalent
 - ✦ 96-well storage plates, 0.8 ml (MIDI plates)

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- ✦ PCR eight-tube strips
- ✦ Solution basin, PVC, RNase/DNase-free (trough)
- ✦ Thermal microplate sealer

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Post-Amp Equipment and Materials

- 1 **96-Well or 384-Well Real-Time PCR Detection System**—One Real-Time PCR Detection system is required. The system must meet the following performance specifications. In addition, the system must have a heated lid and support a 96-well or 384-well format.
 - Minimum Temperature Range: 37°C to 95°C
 - Accuracy/Uniformity: $\pm 0.5^{\circ}\text{C}$
 - Dye: FAM
 - 2 **Thermal Cycler**—One thermal cycler is required. The thermal cycler must have a heated lid and meet the following performance specifications:
 - Minimum Temperature Control Range: 4°C to 99°C
 - Control Accuracy: $\pm 0.3^{\circ}\text{C}$
 - 3 **Optical Compression Pad**—An optical compression pad is recommended to provide a proper thermal seal between the PCR thermal cycler and reaction plates when using adhesive plate seals.
 - 4 **Microplate Shaker**—One microplate shaker is required in the post-amp lab area. The plate shaker must meet the following performance specifications:
 - Maximum Mixing Speed: 3000 rpm
 - Mixing Speed Range: 200 to 3000 rpm
 - 5 **Tabletop Centrifuge**—One tabletop centrifuge is required. (A separate centrifuge is required in the pre-amp area.) Any plate centrifuge that attains the designated speeds of the protocol (280 to $2400 \times g$) is acceptable.
 - 6 **Heat Block**—One heat block for tubes is required. The heat block must meet the following performance specifications.
 - Temperature Range: Ambient $+5^{\circ}\text{C}$ to 99°C
 - Temperature Regulation: $\pm 0.1^{\circ}\text{C}$ at 37°C ; $\pm 0.4^{\circ}\text{C}$ at 60°C
 - 7 **Magnetic Stand**—One magnetic stand for a 96-well plate is required. Better performance is seen when the magnets are on the side of the stand and not on the bottom.
 - 8 **Precision Pipettes**—One set of precision pipettes is required. (A separate set is required in the pre-amp area.) The use of precision pipettes is required to ensure accurate reagent and sample delivery. If they are calibrated regularly and are accurate within 5% of stated volume, single-channel or multi-channel pipettes can be used.
 - 9 **Gel Electrophoresis Supplies**—Gel electrophoresis supplies and apparatus are required along with an appropriate staining method to visualize PCR products in the gel.
 - 10 **Consumables**—The following consumables are required.
 - 96-well skirted PCR plates, 0.2 ml, polypropylene, or equivalent
 - 96-well storage plates, 0.8 ml (MIDI plates)
-  **NOTE**
Make sure that the 96-well plate is fit compatible with the magnetic stand.
- 2-4% TBE Agarose gel

- ▶ 100 bp DNA molecular weight marker
- ▶ DNA loading dye
- ▶ Conical tubes, 15 ml
- ▶ Microcentrifuge tubes
- ▶ PCR eight-tube strips
- ▶ Solution basins, PVC, RNase/DNase-free (trough)
- ▶ Adhesive plate seals
- ▶ Aerosol resistant pipette tips

The PCR process is commonly used in the laboratory to amplify specific DNA sequences. Unless proper laboratory hygiene is used, PCR products can contaminate reagents, instrumentation, and genomic DNA samples, causing inaccurate and unreliable results. PCR product contamination can shut down lab processes and significantly delay normal operations.

Make sure that the lab is set up appropriately to reduce the risk of PCR product contamination:

▶ **Physically Separate Pre- and Post-Amp Areas**

- ▶ Physically separate laboratory space where pre-amp processes are performed (DNA extraction, quantification, and normalization) from the laboratory space where post-amp processes are performed.
- ▶ Never use the same sink to wash pre-amp and post-amp troughs.
- ▶ Never share the same water purification system for pre-amp and post-amp processes.
- ▶ Store all supplies used in the protocols in the pre-amp area, and transfer to the post-amp area as needed.

▶ **Use Dedicated Equipment and Supplies**

- ▶ Dedicate separate full sets of equipment and supplies (pipettes, centrifuges, oven, heat block, etc.) to pre-amp and post-amp lab processes, and never share between processes.
- ▶ Dedicate separate storage areas (freezers and refrigerators) to pre-amp and post-amp consumables.

Because the pre- and post-amplification reagents are shipped together, it is important to unpack the reagents in the pre-amp lab area, and then move the post-amp reagents to the proper post-amp storage area.

To prevent PCR product contamination, it is important to establish lab procedures and follow best practices. Illumina recommends daily and weekly cleaning of lab areas using 0.5% Sodium Hypochlorite (10% Bleach).

CAUTION

To prevent sample or reagent degradation, make sure that all vapors from the cleaning solution have fully dissipated before beginning any processes.

Daily Cleaning of Pre-Amp Area

A daily cleaning of the pre-amp area using a 0.5% Sodium Hypochlorite (10% Bleach) solution helps to eliminate PCR product that has entered the pre-amp area.

Identify pre-amp areas that pose the highest risk of contamination, and clean these areas with a 0.5% Sodium Hypochlorite (10% Bleach) solution before beginning any pre-amp processes. High-risk areas might include, but are not limited to, the following items:

- ▶ Bench tops
- ▶ Door handles
- ▶ Refrigerator/freezer door handles
- ▶ Computer mouse
- ▶ Keyboards

Daily Cleaning of Post-Amp Area

Reducing the amount of PCR product in the post-amp area helps reduce the risk of contamination in the pre-amp area. Daily cleaning of the post-amp area using a 0.5% Sodium Hypochlorite (10% Bleach) solution helps achieve this.

Identify post-amp areas that pose the highest risk of contamination, and clean these areas with a 0.5% Sodium Hypochlorite (10% Bleach) solution daily. High-risk areas might include, but are not limited to, the following items:

- ▶ Thermal cyclers
- ▶ Bench space used to process amplified DNA
- ▶ Door handles
- ▶ Refrigerator/freezer door handles
- ▶ Computer mouse
- ▶ Keyboards

Weekly Cleaning of All Lab Areas

- ▶ Once a week, perform a thorough cleaning of the pre-amp and post-amp areas using 0.5% Sodium Hypochlorite (10% Bleach).
- ▶ Clean all bench tops and laboratory surfaces.
- ▶ Clean all instruments that are not cleaned daily.
- ▶ Thoroughly mop lab floors.
- ▶ Make sure that personnel responsible for weekly cleaning are properly trained on prevention of PCR product contamination.

Items Fallen to the Floor

The floor is contaminated with PCR product transferred on the shoes of individuals coming from the post-amp area; therefore, anything falling to the floor must be treated as contaminated.

- ▶ Disposable items that have fallen to the floor, such as empty tubes, pipette tips, gloves, lab coat hangers, must be discarded.
- ▶ Non-disposable items that have fallen to the floor, such as a pipette or an important sample container, must be immediately and thoroughly cleaned with a 0.5% Sodium Hypochlorite (10% Bleach) solution to remove PCR product contamination.
- ▶ Clean any lab surface that has come in contact with the contaminated item. Individuals handling anything that has fallen to the floor, disposable or non-disposable, must discard their lab gloves and put on a new pair.

Adhere to the following recommendations when preparing libraries for sequencing using this protocol.

Ensuring Consistency

- ▶ **Use multi-channel pipettes**—To ensure consistency across samples, use a multi-channel pipette where possible. Calibrate pipettes regularly.
- ▶ **Consistency for smaller sample preparations**—Each reagent tube supplied with the kit contains sufficient volume to generate results using manual pipettes and reagent troughs following standard laboratory techniques. To ensure accurate reagent volume delivery, single pipette reagent into each well or multi-channel pipette from reagents aliquoted into a PCR 8-tube strip.

Handling Magnetic Beads

- ▶ **Use at room temperature**—Allow the beads to reach room temperature prior to use.
- ▶ **Vortex until well-suspended**—Immediately prior to use, vortex the beads until they are well-suspended and the color appears homogeneous.
- ▶ **Mix samples thoroughly**—After adding the beads to the samples, mix thoroughly by pipetting up and down ten times. Illumina also recommends using a shaker to thoroughly mix samples.
- ▶ **Allow maximum binding**—For best results, incubate the bead/sample mixtures at room temperature for the entire duration indicated in the protocol.
- ▶ **Slowly aspirate cleared solution**—After placing the plate on the magnetic stand, wait for the solution to clear before proceeding. Keep the plate on the magnetic stand when slowly aspirating cleared solution, taking care not to disturb the separated beads.

Avoiding Cross-Contamination

- ▶ **Change tips between dispensing reagents and samples**—Always use fresh pipette tips between dispensing reagents and samples.
- ▶ **Mix plates as directed**—Mix samples with a multi-channel pipette and centrifuge the plate when indicated. Do not vortex the plates.
- ▶ **Use aerosol-resistant tips**—Using aerosol-resistant pipette tips reduces the risk of amplicon carry-over and sample-to-sample cross-contamination.

Washing with 80% Ethanol During the PCR Clean-Up Step

- ▶ **Prepare fresh 80% ethanol**—Always prepare fresh 80% ethanol for wash steps. Ethanol can absorb water from the air and impact results.
- ▶ **Remove all ethanol from wells**—Make sure that all ethanol is removed from the bottom of the wells, as it may contain residual contaminants. Use a P20 multi-channel pipette to remove residual ethanol and accelerate drying.
- ▶ **Allow complete evaporation**—Allow at least five minutes of drying time off the magnetic stand at room temperature for complete evaporation. Residual ethanol can impact the performance of subsequent reactions.

DNA Input Requirements

- ▶ **The Illumina Praxis Extended RAS Panel is designed to be used with genomic DNA** extracted from formalin-fixed paraffin embedded (FFPE) colorectal cancer (CRC) tissue samples.
- ▶ **Assessing DNA Quality**—The Praxis Extended RAS Panel workflow starts with a FFPE DNA qualification step to determine if the DNA is of sufficient quality for use with the

Praxis Extended RAS Panel, by measuring the amplifiable amounts of DNA in the sample relative to a control template. The qualification test determines a difference in quantification cycle (delta Cq value) for each sample relative to the control template.

Genomic DNA (gDNA) samples giving a delta Cq value less than or equal to 5.00 can be used for the Praxis Extended RAS Panel. Illumina does not support the use of gDNA samples giving a delta Cq value of greater than 5.00 with the Praxis Extended RAS Panel.

Quality Controls

Good laboratory practices mandate that the Panel Control sample and a negative (no-template) control sample be included in every run.


Illumina Praxis Extended RAS Panel Acronyms

Acronym	Definition
AMP	AMplification Plate
CLP	CLean-up Plate
DAI	Diluted Amplicon Library
FPU	Filter Plate Unit
HYB	HYBridization Plate
LNP	Library Normalization Plate
NTC	Negative/No Template Control
PAL	Pooled Amplicon Library
SCP	StoraGe Plate

Assay Workflow

The following diagram illustrates the Illumina Praxis Extended RAS Panel workflow. Safe stopping points are marked between steps.

Illumina Praxis Extended RAS Panel Workflow

Qualification of FFPE DNA Samples

During this step, a qPCR reaction will be performed to determine the amplifiability of your FFPE-extracted gDNA samples. By comparing the amplifiability of FFPE DNA relative to that of the Quality Control Template non-FFPE reference gDNA, a ΔC_q value can be calculated for each sample and used to predict its performance in the Praxis Extended RAS Panel assay. The exact amount of FFPE DNA input will vary according to the quality of the extracted DNA.

Estimated Time

- ▶ Total duration: 3 hours
- ▶ Hands-on: 60 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
QC Primers	4 tubes	-25°C to -15°C	Illumina	6
QC Template	1 tube	-25°C to -15°C	Illumina	6
qPCR Master Mix	2 tubes	-25°C to -15°C	Illumina	6
Genomic DNA	As needed	-25°C to -15°C	User	—
RNase/DNase-free Water	As needed	15°C to 30°C	User	—
96 or 384-well Optical PCR Plate and Seal	1 plate	15°C to 30°C	User	—
PCR Eight-tube Strips	As needed	15°C to 30°C	User	—
Microcentrifuge Tubes	As needed	15°C to 30°C	User	—

- 1 Remove the QC Template, QC Primers, qPCR Master Mix, and genomic DNA from -25°C to -15°C storage and thaw at room temperature for up to 30 minutes.
- 2 Place thawed tubes on ice.

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Procedure

- 1 Determine the plate layout of the qPCR reaction.
For 10 samples, Illumina recommends the layout in Figure 2 for a 384-well format. A similar layout can be adapted for a 96-well format. Test QC Template, NTC, and each sample in triplicate.

Suggested qPCR Plate Layout for 10 Samples

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A																								
B																								
C		Sample 1	Sample 1	Sample 1		Sample 7	Sample 7	Sample 7																
D																								
E		Sample 2	Sample 2	Sample 2		Sample 8	Sample 8	Sample 8																
F																								
G		Sample 3	Sample 3	Sample 3		Sample 9	Sample 9	Sample 9																
H																								
I		Sample 4	Sample 4	Sample 4		Sample 10	Sample 10	Sample 10																
J																								
K		Sample 5	Sample 5	Sample 5		QCT	QCT	QCT																
L																								
M		Sample 6	Sample 6	Sample 6		NTC	NTC	NTC																
N																								
O																								
P																								

QCT = QC Template

NTC = No Template Control

- Add 5 µl of QC Template to 495 µl of RNase/DNase-free water in a microcentrifuge tube.
- Vortex the dilution to mix the sample.
- Add 1.5 µl of extracted genomic DNA to 148.5 µl of RNase/DNase-free water in a PCR eight-tube strip to make a 100-fold dilution.

**NOTE**

Illumina recommends performing sample dilutions in a PCR eight-tube strip to facilitate multi-channel pipetting.

- Using a P200, pipette up and down 10 times to mix the dilutions.
- Add 150 µl of the diluted QC Template to an unused well in the PCR eight-tube strip consistent with the plate layout.
- To another unused well, add 150 µl of nuclease-free water for use as NTC sample.
- Prepare the qPCR Master Mix reaction per Table 11 (384-well format) or Table 12 (96-well format). The master mix contains extra volume.

Table 11 qPCR Master Mix Reactions (384-well Format)

Consumable	µl per Well	µl for 36 Reactions
qPCR Master Mix	5.0 µl	205 µl
QC Primers	0.8 µl	33 µl
RNase/DNase-free Water	2.2 µl	90 µl

qPCR Master Mix Reactions (96-well Format)

Consumable	μl per Well	μl for 36 Reactions
qPCR Master Mix	10.0 μl	410 μl
QC Primers	1.6 μl	66 μl
RNase/DNase-free Water	4.4 μl	180 μl

- 9 Mix reaction mix gently but thoroughly.
- 10 Place the reaction mix on ice and protect it from light until use.
- 11 Aliquot master mix into a trough or eight-tube strip to aid in dispensing with multi-channel pipettor.
- 12 Add 8 μl (384-well format) or 16 μl (96-well format) of the master mix to each well of the qPCR plate. Take care to pipette accurately into the wells as small variations will affect the assay.
- 13 Add 2 μl (384-well format) or 4 μl (96-well format) of the QC Template dilution, the sample dilutions, or nuclease-free water to each well of the plate. Take care to pipette accurately into the wells as small variations will affect the assay.
- 14 Seal the plate with the appropriate seal for your instrument, taking care to avoid cross-contamination and to avoid smudging the surface of the lids. Verify that the seal is oriented correctly and uniformly sealed.
- 15 Place the plate on an adapter (if needed) and centrifuge the plate to 1000 × g for 1 minute.
- 16 Ensure that the seal is free of any liquid or dust, place the plate on the qPCR machine in the correct orientation, then close the lid and run the following qPCR thermal profile:
 - ▶ Hot Start
 - ▶ 50°C for 2 minutes
 - ▶ 95°C for 10 minutes
 - ▶ x40
 - ▶ 95°C for 30 seconds
 - ▶ 57°C for 30 seconds
 - ▶ 72°C for 30 seconds
 - ▶ Confirm that the instrument captures images after the 72°C step.

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- 17 After the final step, the thermal cycler analyzes the samples. Ensure that the QC Template produces expected amplification curves (see green amplification curves in Figure 3). On most systems, the QC Template should amplify at a C_q of approximately 15–22 cycles. Exclude replicates from a triplicate group that are > 0.5 C_q different from the rest of the group. At least two of the three replicates less than or equal to 0.5 C_q apart must be included in the final calculation. If these criteria are not met, the qualification process should be repeated for those samples.

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
- 18 Replicates exhibiting abnormal amplification curves should be excluded (see blue amplification curves in Figure 3).

Qualification of Samples* by qPCR

* DNA samples shown in Figure 3 were qualified for use in creating libraries for next-generation sequencing using real-time PCR following parameters stated in the text. The graph shows qPCR amplification curves generated by these samples.

Purple Line—Cq Threshold

Blue Curves—DNA replicates with abnormal amplification profiles that should be excluded from dCq calculations

Green Curves—DNA replicates with expected amplification curves that can be used for delta Cq calculations

- 19 Subtract the average Cq for the QC Template from the average Cq for each sample (Sample Average Cq – QC Template Average Cq) to yield the delta Cq values for each sample. Record the delta Cq values and the sample dilution factors. The sample dilution factors are based on the recommendations in Table 16.
- 20 Make sure that amplification of the NTC occurs at least 10 cycles after QC Template amplification.

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dCq Calculation Examples

Example #	Replicate 1	Replicate 2	Replicate 3	Final values included in Cq calculations	Average Cq	dCq*
1	24.80	25.11	25.06	All	24.99	6.87
Actions:	Does not meet dCq requirement (≤ 5.00). Sample does not qualify for the assay.					
2	17.22	18.58	17.91	None	N/A	N/A
Actions:	Remove replicates >0.5 Cq different from the rest of triplicates (all three). Sample qualification process must be repeated as a minimum of two valid replicates are required.					
3	22.45	22.63	23.52	22.45 & 22.63	22.54	4.42
Actions:	Remove replicate which is >0.5 Cq different from the rest of triplicates. Calculate dCq using two valid replicate values. dCq requirement is met (≤ 5.00).					
4	21.56	21.73	21.48	All	21.59	3.46
Actions:	Sample passed with three replicate values that are not >0.5 Cq apart. dCq requirement is met (≤ 5.00).					
5	18.21	18.50	19.04	All	18.70	0.58
Actions:	Sample passed with three replicate values that are not >0.5 Cq apart. dCq requirement is met (≤ 5.00).					
QC Template (Example 1)	17.86	18.11	18.40	All	18.12	0
Actions:	QC Template passes with three valid replicate values that are not >0.5 Cq apart. QC Template passes QC (average Cq ≥ 15 and ≤ 22).					
	Control considerations					
NTC (negative control)	26.00	25.64	25.83	All	25.82	7.7
Actions:	Entire sample qualification fails. Amplification of the NTC must occur ≥ 10 cycles after QC Template amplification. Repeat qualification process for all samples, QC Template, and NTC.					
QC Template (Example 2)	9.85	9.73	9.87	All	9.82	0
Actions:	QC Template fails. QC Template average Cq must be ≥ 15 and ≤ 22 . Repeat sample qualification process for all samples, QC Template, and NTC.					

*dCq is calculated by [sample Average Cq - QC Template (Example 1) Average Cq]

SAFE STOPPING POINT

Illumina recommends storing extracted genomic DNA between 2°C to 8°C for a maximum of 28 days or storing frozen between -15°C to -25°C for a maximum of 161 days. Genomic DNA stored frozen is stable for a maximum of 2 freeze/thaw cycles.

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After you have populated the samples table, you can export the sample information to an external file, and import the file for another run.

Enter Samples Manually

- 1 Click the plate region A1 or E7 to enter sample information in the samples table.
Region A1 is equivalent to the upper left quadrant of wells in the Hyb plate layout, and includes samples for Oligo Pools A and B.
Region E7 is equivalent to the lower right quadrant of wells in the Hyb plate layout, and includes samples for Oligo Pools A and B.
- 2 In the Pool A table, enter a unique sample name in the Sample Name field for each well you use.
Use alphanumeric characters, dashes, or underscores.
The sample name automatically populates the corresponding well in the Pool B table.
Dual index adapters for each Index Read, specific to the well, display as you enter a sample name.
- 3 [Optional] Although the default layout is recommended, to change the position of control samples, use the Positive Control or Negative Control drop down menu.
- 4 Click the print  icon to display the plate layout.
- 5 Select **Print** to print the plate layout as a reference for preparing libraries.
- 6 Click **Close**.
- 7 Click **Save Run**.

Import Samples

- 1 Click the plate region A1 or E7.
- 2 Click the **Export**  icon for an example template.

- 3 Enter a unique sample name in the template for up to 10 samples and save the file. Use alphanumeric characters, dashes, or underscores.

[Data] Sample Names as a 4 * 3 matrix		
Sample1	Sample5	Sample9
Sample2	Sample6	Sample10
Sample3	Sample7	Positive_Control
Sample4	Sample8	Negative_Control

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Hybridization of Oligonucleotide Pools

During this step, the Extended RAS oligonucleotide pools containing upstream and downstream oligonucleotides specific to the RAS genes are hybridized to genomic DNA samples.

WARNING

This set of reagents contains potentially hazardous chemicals. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Wear protective equipment, including eye protection, gloves, and laboratory coat appropriate for risk of exposure. Handle used reagents as chemical waste and discard in accordance with applicable regional, national, and local laws and regulations. For additional environmental, health, and safety information, see the SDS at support.illumina.com/sds.html.

NOTE

Illumina does not support the use of gDNA samples giving a delta Cq value of greater than 5.00.

Estimated Time

- ▶ Total duration: 3–18 hours (overnight)
- ▶ Hands-on: 15 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
Oligo Pool A - Extended RAS	1 tube	-25°C to -15°C	Illumina	6
Oligo Pool B - Extended RAS	1 tube	-25°C to -15°C	Illumina	6
Hybridization Buffer	1 tube	-25°C to -15°C	Illumina	1A
Genomic DNA	20 µl	-25°C to -15°C	User	—
Panel Control - Extended RAS	1 tube	-25°C to -15°C	Illumina	6
TE Buffer	As needed	15°C to 30°C	User	—
96-well skirted PCR plate	1 plate	15°C to 30°C	User	—
Thermal Microplate Sealer or Adhesive Aluminum Foil	2 seals	15°C to 30°C	User	—
Sterile Troughs	As needed	15°C to 30°C	User	—

Preparation

- 1 Determine genomic DNA samples to be used in the assay. Samples giving a delta Cq value less than or equal to 5.00 can be used for the Praxis Extended RAS Panel.
- 2 [Optional] If extracted DNA samples do not meet the delta Cq threshold, additional tissue sections can be used for DNA extraction if available.
- 3 Remove the Oligo Pool A, Oligo Pool B, Hybridization Buffer, genomic DNA samples, and Panel Control from -25°C to -15°C storage and thaw at room temperature.
- 4 Vortex the Oligo Pool A, Oligo Pool B, and Hybridization Buffer vigorously to make sure that all precipitates have dissolved, then briefly centrifuge the tubes to collect liquid.

**NOTE**

Before using Hybridization Buffer, hold the tube in front of a light and visually inspect to make sure that all precipitates have completely dissolved.

- 5 Set a 96-well heat block to 95°C.
- 6 Pre-heat an incubator to 37°C to prepare for the extension-ligation step.
- 7 Create the run and enter sample plate information in the Praxis Extended RAS Panel analysis module. Verify the location of positive and negative controls and then print sample plate layout. Illumina recommends processing samples in batches no larger than 12 samples/24 libraries.

**NOTE**

Use of controls enables Illumina Technical Support to provide effective troubleshooting assistance. Illumina Technical Support will not provide assistance unless these control reactions were included in the run.

- 1 Label a new 96-well PCR plate "HYB_Plate_ID".
See Figure 4 for the suggested HYB plate layout for the following steps.

HYB Plate Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	Sample 5	Sample 9	Sample 1	Sample 5	Sample 9						
B	Sample 2	Sample 6	Sample 10	Sample 2	Sample 6	Sample 10						
C	Sample 3	Sample 7	Positive Control	Sample 3	Sample 7	Positive Control						
D	Sample 4	Sample 8	Negative Control	Sample 4	Sample 8	Negative Control						
E							Sample 11	Sample 15	Sample 19	Sample 11	Sample 15	Sample 19
F							Sample 12	Sample 16	Sample 20	Sample 12	Sample 16	Sample 20
G							Sample 13	Sample 17	Positive Control	Sample 13	Sample 17	Positive Control
H							Sample 14	Sample 18	Negative Control	Sample 14	Sample 18	Negative Control

Extended RAS-Oligo Pool A

Extended RAS-Oligo Pool B

- 2 Using the plate layout in Figure 4, add 10 µl of TE Buffer to the wells indicated as Negative Control (No Template Control) for Oligo Pool A and Oligo Pool B.
- 3 For each 12 sample/ 24 library preparation, add 2.5 µl of Panel Control to 22.5 µl of TE Buffer in a microcentrifuge tube. Mix the diluted Panel Control thoroughly.
- 4 Using the plate layout in Figure 4, add 10 µl of diluted Panel Control prepared in the previous step to each of the two indicated Positive Control wells for Oligo Pool A and Oligo Pool B.
- 5 Determine if sample dilution is required based on Table 16. If dilution is necessary, take 15 µl of genomic DNA extracted from FFPE samples and dilute in TE buffer according to Table

16. For example, for a sample with delta C_q of -1.50, add 15 µl sample to 15 µl of TE Buffer. For a sample with delta C_q of -2.50, add 15 µl sample to 45 µl of TE Buffer.

Delta C_q Dilution Table

Delta C _q	< -1.51	-1.50 to -0.51	-0.50 to 0.50	0.51 to 1.50	1.51 to 5.00	> 5.00
Dilution	4x	2x	No dilution	No dilution	No dilution	Not supported

- 6 Using the plate layout in Figure 4, add 10 µl of each genomic DNA sample as prepared in step 5 to the Oligo Pool A wells of the **HYB** plate, starting with columns 1 or 7. Repeat this process for the Oligo Pool B wells of the **HYB** plate, starting with columns 4 or 10.

NOTE

Verify that DNA sample layout and the positions of positive and negative controls match the plate layout in Figure 4.

- 7 Using a multi-channel pipette, add 5 µl of Oligo Pool A to all wells containing genomic DNA, Panel Control, and Negative Control indicated for Oligo Pool A in columns 1-3 or 7-9 in the **HYB** plate. Then add 5 µl of Oligo Pool B to all wells containing genomic DNA, Panel Control, and Negative Control indicated for Oligo Pool B in columns 4-6 or 10-12 in the **HYB** plate. Change tips after each column to avoid cross-contamination.
- 8 Using a multi-channel pipette, add 40 µl of Hybridization Buffer to each sample in the **HYB** plate. Gently pipette up and down 3–5 times to mix. Change tips after each column to avoid cross-contamination.

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NOTE

Do not mix the Oligo Pools and Hybridization Buffer for storage. If combined, the Oligo Pools become unstable, even when stored frozen.

- 9 Seal the **HYB** plate with a thermal microplate sealer. If a thermal sealer is not available, use an adhesive aluminum foil seal and secure the seal with a sealing roller or wedge, especially on the perimeter of the plate.
- 10 Centrifuge at 1000 × g at room temperature for 1 minute.
- 11 Place the **HYB** plate in the pre-heated block at 95°C and incubate for 1 minute.
- 12 Change the temperature of the same heat block to 40°C and incubate for 3–18 hours.

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This process removes unbound oligonucleotides from genomic DNA using a filter capable of size selection. Two wash steps using Stringent Wash Buffer ensure complete removal of unbound oligonucleotides. A third wash step using Universal Wash Buffer removes residual Stringent Wash Buffer and prepares samples for the extension-ligation step.

WARNING

This set of reagents contains potentially hazardous chemicals. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Wear protective equipment, including eye protection, gloves, and laboratory coat appropriate for risk of exposure. Handle used reagents as chemical waste and discard in accordance with applicable regional, national, and local laws and regulations. For additional environmental, health, and safety information, see the SDS at support.illumina.com/sds.html.


Estimated Time

- ▶ Total duration: 25 minutes
- ▶ Hands-on: 25 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
Extension-Ligation Mix	1 tube	-25°C to -15°C	Illumina	1A
Stringent Wash Buffer	1 bottle	2°C to 8°C	Illumina	3A
Universal Wash Buffer	1 tube	2°C to 8°C	Illumina	3A
TruSeq™ Custom Amplicon Dx Filter Plate	1 plate	15°C to 30°C	Illumina	5
Adapter collar	1 plate	15°C to 30°C	User	—
MIDI plate	1 plate	15°C to 30°C	User	—
Troughs	As needed	15°C to 30°C	User	—

- 1 Remove Extension-Ligation Mix from -25°C to -15°C storage and thaw at room temperature.
- 2 Extension-Ligation Mix is used in the Extension-Ligation step and takes approximately 20 minutes to thaw.
- 3 Remove Stringent Wash Buffer and Universal Wash Buffer from 2°C to 8°C storage and set aside at room temperature.
- 4 Assemble the filter plate assembly unit (FPU) in the following order (from bottom to top):

- A Lid
- B Filter plate
- C Adapter collar
- D MIDI plate

- 5 Label the filter plate **FPU_Plate_ID**.
- 6 Pre-wash the filter plate membrane as follows:
 - a Using a multi-channel pipette, add 50 μ l of Stringent Wash Buffer to each well.
 - b Cover the **FPU** plate with the filter plate lid and keep it covered during each centrifugation step.
 - c Centrifuge the **FPU** at $2400 \times g$ at room temperature for 5 minutes.
 - d If the wash buffer does not drain completely, centrifuge the filter plate again at $2400 \times g$ at room temperature for 5 minutes.

NOTE

Illumina strongly recommends keeping spare filter plates (FC-130-1006) on hand as general lab supplies.

Procedure

- 1 After the hybridization incubation, confirm the heat block has cooled to 40°C.

NOTE
If the heat block fails to cool to 40°C overnight, library preparation will need to be repeated.
- 2 Remove the **HYB** plate from the heat block and centrifuge at $1000 \times g$ at room temperature for 1 minute to collect condensation.

- 3 Using a multi-channel pipette set to 65 μ l, transfer the entire volume of each sample onto the center of the corresponding pre-washed wells of the filter plate. Change tips after each column to avoid cross-contamination.

- 4 Cover the filter plate with the lid and centrifuge at $2400 \times g$ at room temperature for 5 minutes.

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- 5 Wash the filter plate as follows:

- a Using a multi-channel pipette, add 50 μ l of Stringent Wash Buffer to each sample well. Changing tips between columns is not required if you use care to avoid cross-contamination.
- b Cover the filter plate with the lid and centrifuge at $2400 \times g$ at room temperature for 5 minutes.
- c If the wash buffer does not drain completely, centrifuge the filter plate again at $2400 \times g$ at room temperature for 5 minutes.

- 6 Repeat the wash as follows:

- a Using a multi-channel pipette, add 50 μ l of Stringent Wash Buffer to each sample well. Changing tips between columns is not required if you use care to avoid cross-contamination.
- a Cover the filter plate with the lid and centrifuge at $2400 \times g$ at room temperature for 5 minutes.
- b If the wash buffer does not drain completely, centrifuge the filter plate again at $2400 \times g$ at room temperature for 5 minutes.

- 7 Discard all the flow-through (containing formamide) collected up to this point in an appropriate hazardous waste container, then reassemble the FPU. The same MIDI plate can be re-used for the rest of the pre-amplification process.

- 8 Using a multi-channel pipette, add 45 μ l of Universal Wash Buffer to each sample well. Changing tips between columns is not required if you use care to avoid cross-contamination.

- 9 Cover the filter plate with the lid and centrifuge at $2400 \times g$ at room temperature for 5 minutes.

NOTE

Make sure that all liquid has drained after centrifugation. Repeat centrifugation if necessary. Residual wash buffer may inhibit subsequent enzymatic reactions.

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This process connects the hybridized upstream and downstream oligonucleotides. A DNA polymerase extends from the upstream oligonucleotide through the targeted region, followed by ligation to the 5' end of the downstream oligonucleotide using a DNA ligase. This results in the formation of products containing the targeted regions of interest flanked by sequences required for amplification.

Estimated Time

- ▶ Total duration: 50 minutes
- ▶ Hands-on: 5 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
Extension-Ligation Mix	1 tube	-25°C to -15°C	Illumina	1A
Adhesive aluminum foil seal	1 seal	15°C to 30°C	User	—
Troughs	As needed	15°C to 30°C	User	—

- 1 Using a multi-channel pipette, add 45 µl of Extension-Ligation Mix to each sample well of the filter plate. The Extension-Ligation reaction takes place on the filter plate membrane. Changing tips between columns is not required if you use care to avoid cross-contamination.
- 2 Seal the filter plate with adhesive aluminum foil, and then cover with the lid to secure the foil during incubation.
- 3 Incubate the entire FPU assembly in the pre-heated 37°C incubator for 45 minutes.
- 4 While the FPU plate is incubating, prepare the AMP (Amplification Plate) as described in the following section.

PCR Amplification

In this step, the extension-ligation products are amplified using primers that add index sequences for sample multiplexing, as well as common adapters required for cluster generation.

Estimated Time


- ✦ Total duration: ~105 minutes
- ✦ Hands-on: 30 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
PCR Master Mix	1 tube	-25°C to -15°C	Illumina	1A
Index Primers A-H (A501-A508)	1 tube per primer	-25°C to -15°C	Illumina	1A
Index Primers 1-12 (A701-A712)	1 tube per primer	-25°C to -15°C	Illumina	1A
PCR Polymerase	1 tube	-25°C to -15°C	Illumina	1A
Appropriate PCR plate seal	1	15°C to 30°C	User	—
0.05 N NaOH, freshly prepared	5 mL	15°C to 30°C	User	—
96-well skirted PCR plate	1 plate	15°C to 30°C	User	—
Troughs	As needed	15°C to 30°C	User	—

- 1 Prepare fresh 0.05 N NaOH by adding 25 µl of 10 N NaOH to 4975 µl of sterile water.
- 2 Determine the index primers to be used according to the plate graphic printout from Praxis Extended RAS Panel analysis module.
- 3 Remove PCR Master Mix and the appropriate index primers from -25°C to -15°C storage and thaw on a bench at room temperature. Allow approximately 20 minutes to thaw the reagents.
- 4 After the index primers are completely thawed, vortex each tube to mix and briefly centrifuge the tubes in a microcentrifuge. Use 1.7 ml microcentrifuge tubes as adapters for the microcentrifuge.
- 5 Arrange the primers in a rack using the following arrangements:
 - a Arrange Index Primers A-H (A501-A508) primer tubes (white caps, clear solution) vertically, aligned with rows A through H.
 - b Arrange Index Primers 1-12 (A701-A712) primer tubes (orange caps, yellow solution) horizontally, aligned with columns 1 through 12.

Figure 6 Index Plate Fixture

- A Index Primers A-H (A501-A508) (white caps)
- B Index Primers 1-12 (A701-A712) (orange caps)
- C AMP plate

- 6 Label a new 96-well PCR plate "AMP" (Amplification Plate).
- 7 Add index primers to the AMP plate as follows:
 - a Using a multi-channel pipette, add 9 μ l of the selected index primers [A-H (A501-A508)] (clear solution) to the appropriate well in a column of the AMP plate.
 - b Using a multi-channel pipette, add 9 μ l of the selected index primers [1-12 (A701-A712)] (yellow solution) to the appropriate row of the AMP plate. *Tips must be changed after each row to avoid index cross-contamination.*
 - c Remove all the index primer tubes from the working area.

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- 8 Prepare the PCR Master Mix/PCR Polymerase PCR working solution as follows:
 - a For 12 samples/24 libraries, add 21 μ l of PCR Polymerase to 1050 μ l of PCR Master Mix.
 - b Invert the prepared PCR working solution 20 times to mix.

You will add this working solution to the AMP plate in the next section. The PCR working solution is stable at room temperature for 10 minutes and can also be stored on ice if needed during preparation.

NOTE

Always add PCR Polymerase to PCR Master Mix just prior to use. Never store the combined PCR working solution.

Procedure

- 1 When the 45-minute extension-ligation reaction is complete, remove the **FPU** from the incubator. Remove the aluminum foil seal and replace with the filter plate lid. Removing the aluminum foil seal before centrifugation ensures the reaction supernatant will drain into the waste plate effectively.
 - 2 Centrifuge the **FPU** at $2400 \times g$ at room temperature for 2 minutes.
 - 3 Using a multi-channel pipette, add 25 μl of 0.05 N NaOH to each sample well on the filter plate. Ensuring that pipette tips come in contact with the membrane, pipette the NaOH up and down 5–6 times. Tips must be changed after each column.
 - 4 Cover and incubate the filter plate at room temperature for 5 minutes.
 - 5 While the filter plate is incubating, use a multi-channel pipette to transfer 22 μl of the PCR working solution to each well of the AMP plate containing index primers. Change tips between samples.
 - 6 Transfer samples eluted from the filter to the AMP plate as follows:
 - a Set a multi-channel P20 pipette to 20 μl .
 - b Pipette the samples in the first column of the filter plate up and down 5–6 times.
- 

NOTE
Slightly tilt the **FPU** plate to ensure complete aspiration and to avoid air bubbles.
- c Transfer 20 μl from the filter plate to the corresponding column of the **AMP** plate.
 - d Gently pipette up and down 5–6 times to thoroughly combine the DNA with the PCR working solution.
 - e Transfer the remaining columns from the filter plate to the AMP plate in a similar manner. *Tips must be changed after each column to avoid index and sample cross-contamination.*
 - f After all the samples have been transferred, the waste collection MIDI plate of the **FPU** can be discarded. The metal adapter collar should be cleaned and put away for future use.
- 7 Cover the **AMP** plate with the appropriate plate seal and secure with a rubber roller.
 - 8 Centrifuge at $1000 \times g$ at room temperature for 1 minute.
 - 9 Transfer the AMP plate to the post-amplification area.
 - 10 Perform PCR using the following program on a thermal cycler:
 - ▶ 95°C for 3 minutes
 - ▶ 32 cycles of:
 - ▶ 95°C for 30 seconds
 - ▶ 66°C for 30 seconds
 - ▶ 72°C for 60 seconds
 - ▶ 72°C for 5 minutes
 - ▶ Hold at 10°C

NOTE

Illumina strongly recommends using an Optical Compression Pad to provide a proper thermal seal between the PCR thermal cycler and reaction plates and prevent sample evaporation.

NOTE

On the PCR machine, set the reaction volume to 60 μ l and the temperature ramp speed to maximum.

SAFE STOPPING POINT

If not proceeding immediately to PCR Clean-Up, the AMP plate can remain on the thermal cycler overnight, or can be stored at 2°C to 8°C up to 48 hours.

Venty Library Preparation

After PCR:

- 1 Combine 5 μ l of amplified product with 15 μ l of RNase/DNase-free water and DNA loading dye if necessary.
- 2 Run on a 2-4% TBE agarose gel with 100 bp ladder to confirm the presence and brightness of the 250-280 bp library product.
 - Samples that show amplification in both Oligo Pool A and B are considered valid and can be processed through the remainder of the workflow (Figure 7).
 - Samples showing little to no amplification in one Oligo Pool (Figure 8, left) or both Oligo Pools (Figure 8, right) are considered invalid and should not be processed through the remainder of the workflow. If an invalid gel result is observed, library preparation for that sample or samples will need to be repeated.
 - The Negative Control sample should show an Invalid gel result. If the Negative Control sample shows amplification in both Oligo Pool A and B, this indicates potential processing errors, i.e. contamination, pipetting errors, or sample mix-ups. The operator should evaluate the root cause and repeat library preparation as appropriate. If contamination is suspected, library preparation for all samples should be repeated. If pipetting errors or mix-ups are suspected, repeat library preparation for the affected samples.
 - The presence of a 500 bp band is normal and does not affect assay performance.

Figure 7 Valid Gel Results, Example




Figure 8 Invalid Gel Results, Examples

- A 100 bp Ladder
- B Oligo Pool A
- C Oligo Pool B

Figure 9 Sample Outcomes, Example

- A 100 bp Ladder
- B Oligo Pool A
- C Oligo Pool B

PCR Clean-Up

This process uses PCR Clean-Up Beads to purify the PCR products from the other reaction components.

Estimated Time

- ▶ Total duration: 50 minutes
- ▶ Hands-on: 20 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
Elution Buffer	1 tube	15°C to 30°C	Illumina	5
PCR Clean-Up Beads	1.8 ml per 12 samples/24 libraries	2°C to 8°C	Illumina	3B
80% Ethanol, freshly-prepared	15 ml per 12 samples/24 libraries	15°C to 30°C	User	—
96-well MIDI plates	2	15°C to 30°C	User	—
Adhesive plate seal	As needed	15°C to 30°C	User	—
Troughs	As needed	15°C to 30°C	User	—

Precautions

NOTE

Please review the **Precautions** section at the beginning of this protocol regarding the handling of magnetic beads and washing with 80% ethanol during the PCR clean-up.

- 1 Bring the PCR Clean-Up Beads to room temperature.
- 2 Prepare fresh 80% ethanol from absolute ethanol. For 12 samples/24 libraries, prepare 80% ethanol using 12 mL of absolute ethanol and 3 mL of DNase/RNase-free water. Mix thoroughly.

NOTE

Always prepare fresh 80% ethanol for wash steps. Ethanol can absorb water from the air impacting the results.

Procedure

- 1 Centrifuge the AMP plate at 1000 × g at room temperature for 1 minute to collect condensation.
- 2 Label a new MIDI plate "CLP_Plate_ID" (Clean-up Plate).
- 3 Invert PCR Clean-Up Beads 10 times. Vortex vigorously and then invert again 10 times.
- 4 Visually inspect solution to ensure that beads are resuspended well.
- 5 Using a multi-channel pipette, add 55 µl of PCR Clean-Up Beads to each well of the CLP plate.

- 6 Using a multi-channel pipette set to 60 μ l, transfer the entire PCR product from the **AMP** plate to the **CLP** plate. Change tips between samples.
- 7 Seal the **CLP** plate with an adhesive plate seal.
- 8 Shake the **CLP** plate on a microplate shaker at 1800 rpm for 2 minutes.
- 9 Incubate at room temperature (15°C to 30°C) without shaking for 10 minutes.
- 10 Place the plate on a magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
- 11 With the **CLP** plate on the magnetic stand and a multi-channel pipette set to 100 μ l, carefully remove and discard the supernatant. Change tips between samples.

NOTE

If any beads are inadvertently aspirated into the tips, dispense the beads back to the plate and let the plate rest on the magnet for 2 minutes and confirm that the supernatant has cleared.

NOTE

Delays during this step may lead to bead clumping following removal of the supernatant. Proceed immediately to the next step as soon as all supernatant is removed.

- 12 With the **CLP** plate on the magnetic stand, wash the beads with freshly prepared 80% ethanol as follows:
 - a Using a multi-channel pipette, add 200 μ l of freshly prepared 80% ethanol to each sample well. Changing tips is not required if you use care to avoid cross-contamination. You do not need to resuspend the beads at this time.
 - b Incubate the plate on the magnetic stand for a minimum of 30 seconds or until the supernatant is clear.
 - c Carefully remove and discard the supernatant.
 - 13 With the **CLP** plate on the magnetic stand, perform a second ethanol wash as follows:
 - a Using a multi-channel pipette, add 200 μ l of freshly prepared 80% ethanol to each sample well.
 - b Incubate the plate on the magnetic stand for a minimum of 30 seconds or until the supernatant appears clear.
 - c Carefully remove and discard the supernatant.
 - 14 Use a P20 multi-channel pipette set to 20 μ l to remove excess ethanol. Do not disturb or touch the beads.
 - 15 Remove the **CLP** plate from the magnetic stand and allow the beads to air-dry for 5 minutes.
 - 16 Using a multi-channel pipette, add 25 μ l of Elution Buffer to each sample. Changing tips is not required if you use care to avoid cross-contamination.
 - 17 Seal the plate with an adhesive plate seal.
 - 18 Shake the **CLP** plate on a microplate shaker at 1800 rpm for 5 minutes.
- 

NOTE
Ensure all samples are completely resuspended. If there are samples in which the beads are not completely resuspended, gently pipette up and down to resuspend the beads and repeat the previous two steps.
- 19 Place the **CLP** plate on the magnetic stand for a minimum of 2 minutes or until the supernatant is clear.

- 20 Label a new MIDI plate "LNP_Plate_ID" (Library Normalization Plate).
- 21 Using a P20 multi-channel pipette and fine tips, carefully transfer 20 µl of the supernatant from the CLP plate to the LNP plate. Change tips between samples to avoid cross-contamination.

 **NOTE**
If any beads are inadvertently aspirated into the tips, dispense the beads back to the plate and let the plate rest on the magnet for 2 minutes and confirm that the supernatant has cleared.
- 22 [Optional] Transfer the remaining 5 µl of supernatant from the CLP plate to a new plate and label the plate with a run name and date. Store this plate at -25°C to -15°C until completion of the sequencing run and data analysis. The cleaned up PCR products can be used for troubleshooting efforts in the event of sample failures.
- 23 Seal the LNP plate with an adhesive plate seal, and then centrifuge at 1000 × g at room temperature for 1 minute to ensure all the supernatant is at the bottom of the well.

SAFE STOPPING POINT

After PCR Clean-Up, the plate is stable for up to 3 hours at 2°C to 8°C.

Library Normalization

This process normalizes the quantity of each library to ensure equal library representation in the pooled sample.

Estimated Time

- ▶ Total duration: 1 hour 20 minutes
- ▶ Hands-on: 30 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
Library Normalization Diluent	1 tube	-25°C to -15°C	Illumina	1B
Library Beads	1 tube	2°C to 8°C	Illumina	3B
Library Normalization Wash	2 tubes	2°C to 8°C	Illumina	3B
Library Storage Buffer	1 tube	15°C to 30°C	Illumina	5
0.1 N NaOH, freshly prepared	3 ml	15°C to 30°C	User	—
96-well skirted PCR plate	1 plate	15°C to 30°C	User	—
15 ml conical tube	1 tube	15°C to 30°C	User	—
Adhesive plate seal	As needed	15°C to 30°C	User	—

WARNING

This set of reagents contains potentially hazardous chemicals. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Wear protective equipment, including eye protection, gloves, and laboratory coat appropriate for risk of exposure. Handle used reagents as chemical waste and discard in accordance with applicable regional, national, and local laws and regulations. For additional environmental, health, and safety information, see the SDS at support.illumina.com/sds.html.

Preparation

- 1 Prepare fresh 0.1 N NaOH by adding 30 µl of 10 N NaOH to 2970 µl of RNase/DNase-free water.
- 2 Remove Library Normalization Diluent from -25°C to -15°C storage and bring to room temperature. Use a 20°C to 25°C water bath as needed.

NOTE

Library Normalization Diluent might form visible precipitates or crystals. Before use, vortex vigorously, and then hold the tube in front of a light and visually inspect to make sure that all precipitate has completely dissolved.

- 3 Remove Library Beads and Library Normalization Wash from 2°C to 8°C storage and bring to room temperature.
Use a 20°C to 25°C water bath as needed.

- 4 Vortex Library Beads vigorously for 1 minute with intermittent inversion until the beads are resuspended and no pellet is found at the bottom of the tube when the tube is inverted.

Procedure

- 1 For 12 samples/24 libraries, add 1575 μ l of Library Normalization Diluent to a fresh 15 ml conical tube.
- 2 To homogenize beads, use a P1000 pipette set to 1000 μ l to resuspend Library Beads thoroughly by pipetting up and down 10 times.

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- 3 For 12 samples/24 libraries, pipette 300 μ l of Library Beads to the tube containing Library Normalization Diluent. Mix well by inverting the tube 15–20 times.

NOTE

A P1000 set to 1000 μ l is required to resuspend the beads completely in step 2. Mix only the specified amounts of Library Normalization Diluent and Library Beads. You must store the remaining Library Normalization Diluent and Library Beads separately at their respective recommended temperatures. To preserve stability, Library Beads should never be frozen or mixed with Library Normalization Diluent if not used immediately.

- 4 Using a multi-channel pipette, add 45 μ l of the combined Library Normalization Diluent/Library Beads working solution to each well of the LNP plate containing libraries. Changing tips between columns is not required if you use care to avoid cross-contamination.
- 5 Seal the LNP plate with an adhesive plate seal.
- 6 Shake the LNP plate on a microplate shaker at 1,800 rpm for 30 minutes.

NOTE

This 30-minute incubation is critical for proper library normalization. Incubations of greater or less than 30 minutes may affect library representation and cluster density.

- 7 Place the LNP plate on a magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
- 8 With the LNP plate on the magnetic stand, using a multi-channel pipette set to 80 μ l, carefully remove and discard the supernatant in an appropriate hazardous waste container.

NOTE

If any beads are inadvertently aspirated into the tips, dispense the beads back to the plate and let the plate rest for 2 minutes or until the supernatant has cleared.

- 9 Remove the LNP plate from the magnetic stand and wash the beads with Library Normalization Wash as follows:
 - a Using a multi-channel pipette, add 45 μ l of Library Normalization Wash to each sample well.
Changing tips between columns is not required if you use care to avoid cross-contamination.
 - b Seal the LNP plate with an adhesive plate seal.
 - c Shake the LNP plate on a microplate shaker at 1,800 rpm for 5 minutes.

- d Place the **LNP** plate on the magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
 - e Carefully remove and discard the supernatant in an appropriate hazardous waste container.
 - 10 Remove the **LNP** plate from the magnetic stand and repeat the wash with Library Normalization Wash as follows:
 - a Using a multi-channel pipette, add 45 μ l of Library Normalization Wash to each well. Changing tips between columns is not required if you use care to avoid cross-contamination.
 - b Seal the **LNP** plate with an adhesive plate seal.
 - c Shake the **LNP** plate on a microplate shaker at 1,800 rpm for 5 minutes.
 - d Place the **LNP** plate on the magnetic stand for a minimum of 2 minutes.
 - e Carefully remove and discard the supernatant in an appropriate hazardous waste container.
 - f Seal the **LNP** plate with an adhesive plate seal.
 - 11 Centrifuge the **LNP** plate at $1000 \times g$ at room temperature for 30 seconds to collect residual wash buffer.
 - 12 Place the **LNP** plate on the magnetic stand.
 - 13 Use a P20 multi-channel pipette set to 20 μ l to remove excess Library Normalization Wash.
 - 14 Remove the **LNP** plate from the magnetic stand and add 30 μ l of 0.1 N NaOH to each well to elute the sample.
 - 15 Seal the **LNP** plate with an adhesive plate seal.
 - 16 Shake the **LNP** plate on a microplate shaker at 1,800 rpm for 5 minutes.
 - 17 During the 5-minute shake, label a new 96-well PCR plate **SGP_Plate_ID**.
 - 18 Add 30 μ l Library Storage Buffer to each well to be used in the **SGP** plate.
 - 19 After the 5-minute shake, ensure that all samples in the **LNP** plate are resuspended. If the samples are not resuspended, gently pipette those samples up and down or lightly tap the plate on the bench to resuspend the beads, then shake for another 5 minutes.
 - 20 [Optional] Centrifuge the **LNP** plate at $1000 \times g$ at room temperature for 30 seconds to collect contents to the bottom of the wells.
 - 21 Place the **LNP** plate on the magnetic stand for a minimum of 2 minutes.
 - 22 Using a multi-channel pipette set to 30 μ l, transfer the supernatant from the **LNP** plate to the **SGP** plate. Gently pipette up and down 5 times to mix.
-  **NOTE**
If any beads are inadvertently aspirated into the tips, dispense the beads back to the plate and let the plate rest for 2 minutes or until the supernatant has cleared.
- 23 Seal the **SGP** plate with an adhesive plate seal and then centrifuge at $1,000 \times g$ at room temperature for 1 minute. Proceed immediately to Library Pooling.

Library Pooling

In preparation for cluster generation and sequencing, equal volumes of normalized library are combined, diluted in hybridization buffer, and heat denatured prior to sequencing on the MiSeqDx. PhiX is used as an internal control for sequencing.

Estimated Time

- ▶ Total duration: 10 minutes
- ▶ Hands-on: 10 minutes

Consumables

Item	Quantity	Storage	Supplied By	Kit Box
Library Dilution Buffer	2 tubes	-25°C to -15°C	Illumina	1B
PhiX Internal Control	2 µl	-25°C to -15°C	Illumina	1B
0.1 N NaOH, freshly prepared	2.5 ml	15°C to 30°C	User	—
1X TE buffer	8 µl	15°C to 30°C	User	—
Microcentrifuge tubes	As needed	15°C to 30°C	User	—
PCR eight-tube strip	1	15°C to 30°C	User	—
Ice bucket	1	15°C to 30°C	User	—


Prepare for Library Pooling

- 1 Set a heat block suitable for 1.5 ml centrifuge tubes to 96°C.
- 2 In an ice bucket, prepare an ice-water bath. Chill the Library Dilution Buffer in the ice-water bath.
- 3 Begin thawing the MiSeq™ Dx SBS Cartridge v3.

Prepare for Reagent Cartridge

The following instructions describe how to thaw the reagent cartridge using a room temperature water bath. This method requires approximately one hour.

- 1 Remove the reagent cartridge from -25°C to -15°C storage.
- 2 Place the reagent cartridge in a water bath containing enough room temperature deionized water to submerge the base of the reagent cartridge up to the water line printed on the reagent cartridge. Do not allow the water to exceed the maximum water line.

- 3 Allow the reagent cartridge to thaw in the room temperature water bath for approximately 60-90 minutes or until completely thawed.
- 4 Remove the cartridge from the water bath and gently tap it on the bench to dislodge water from the base of the cartridge. Dry the base of the cartridge. Make sure that no water has splashed on the top of the reagent cartridge.

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- ▶ DNase, RNase-free water (2475 μ l)
- ▶ Stock 10 N NaOH (25 μ l)

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- ▶ PhiX Internal Control (2 μ l)
- ▶ 1X TE Buffer (8 μ l)

- 4 Combine the following volumes of Diluted PhiX Internal Control and 0.1 N NaOH in a microcentrifuge tube to make Denatured Diluted PhiX Internal Control:

- ▶ Diluted PhiX Internal Control (10 μ l)
- ▶ 0.1 N NaOH (10 μ l)

- 5 Vortex briefly to mix the Denatured Diluted PhiX Internal Control.
- 6 Briefly centrifuge the Denatured Diluted PhiX Internal Control to collect contents.
- 7 Incubate for 5 minutes at room temperature to denature the PhiX Internal Control into single strands.
- 8 Combine the following volume of pre-chilled Library Dilution Buffer and Denatured Diluted PhiX Internal Control to make PhiX/Library Dilution Buffer solution.

- ▶ Denatured Diluted PhiX Internal Control (20 μ l)
- ▶ Pre-chilled Library Dilution Buffer (980 μ l)

NOTE

You can store the the PhiX/Library Dilution Buffer solution up to 3 weeks at -25°C to -15°C as single-use aliquots. After 3 weeks, cluster numbers tend to decrease.

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- 2 Inspect reagents in positions 1, 2, and 4 to make sure that they are fully mixed and free of precipitates.
- 3 Gently tap the cartridge on the bench to reduce air bubbles in the reagents.

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Preparing Samples for Sequencing

- 1 Bring Library Dilution Buffer to room temperature. Vortex Library Dilution Buffer and ensure that all precipitates have dissolved completely.
- 2 Centrifuge the **SGP** plate at 1000 × g at room temperature for 1 minute to collect condensation.
- 3 Label a fresh microcentrifuge tube "**PAL_Plate_ID**" (Pooled Amplicon Library).
- 4 Determine the samples to be pooled for sequencing. A maximum of 12 samples/24 libraries can be pooled for each sequencing run.
- 5 Transfer 10 µl of each library to be sequenced from the **SGP** plate, column by column, to a PCR eight-tube strip. Seal **SGP** with an adhesive plate seal and set aside.
- 6 Combine and transfer the contents of the PCR eight-tube strip into the **PAL** tube. Mix the **PAL** tube thoroughly.
- 7 Label 3–6 fresh microcentrifuge tubes "**DAL_Plate_ID**" (Diluted Amplicon Library).

**NOTE**

Illumina recommends creating 3 **DAL** tubes per **PAL** for future use.

- 8 Add 784 µl of Library Dilution Buffer to each **DAL** tube.
- 9 Add 8 µl of PhiX/Library Dilution Buffer solution to each **DAL** tube. Using the same tip, pipette up and down 3–5 times to rinse the tip and ensure complete transfer.
- 10 Transfer 8 µl of **PAL** to each **DAL** tube containing Library Dilution Buffer. Using the same tip, pipette up and down 3–5 times to rinse the tip and ensure complete transfer.
- 11 Mix the **DAL** tubes by vortexing the tube at top speed.

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- 12 Centrifuge the one **DAL** tube at room temperature for 1 minute to collect contents.
- 13 Using a heat block, incubate the **DAL** tube at 96°C for 2 minutes.

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- 14 After the incubation, invert the **DAL** tube 1–2 times to mix, then immediately place in the ice-water bath.

- 15 Keep the DAL tube in the ice-water bath for 5 minutes.

NOTE

You must perform the heat denaturation step immediately before loading DAL into the MiSeq™ SBS cartridge v3 to ensure efficient template loading on the MiSeq™ Dx flow cell. Replicate DAL tubes must be heat denatured immediately before MiSeqDx loading.

Library Sequencing

The flow cell is washed, dried, and loaded into MiSeqDx, samples are loaded into the reagent cartridge, the reagent cartridge is loaded into MiSeqDx, and the sequencing run is started. The MiSeqDx performs cluster generation, sequencing by synthesis, and data analysis.

Estimated Time

- ▶ Total sequencing duration: ~28 hours
- ▶ Hands-on: ~5 minutes


Consumables

Item	Quantity	Storage	Supplied By	Kit Box
MiSeq™ Dx SBS Cartridge v3	1 pre-filled cartridge	-25°C to -15°C	Illumina	2
MiSeq™ Dx SBS Solution (PR2)	1 bottle	2°C to 8°C	Illumina	4
MiSeq™ Dx Flow Cell	1 flow cell	2°C to 8°C	Illumina	3B

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- 1 Use a separate, clean, and empty 1 ml pipette tip to pierce the foil seal over the reservoir on the MiSeq™ Dx SBS Cartridge v3 labeled **Load Samples**.
- 2 Pipette 600 µl of the sample DAL libraries into the **Load Samples** reservoir. Take care to avoid touching the foil seal while dispensing the sample.
- 3 Check for air bubbles in the reservoir after loading sample. If air bubbles are present, gently tap the cartridge on the bench to release the bubbles.

Load Libraries

- 4 Log in to the MiSeq Operating Software (MOS) using your Local Run Manager user name and password.

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- 5 From the Home screen, select **Sequence**.
- 6 Select a run from the list.
- 7 Confirm the run settings.
- 8 Clean the flow cell.
- 9 Load the flow cell.
- 10 Empty the waste bottle and load the MiSeq[™] Dx SBS Solution (PR2) bottle.
- 11 Load the reagent cartridge.
- 12 The run automatically starts after passing pre-run checks.
- 13 Perform a post-run wash.

[This page intentionally left blank]

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: www.illumina.com
Email: techsupport@illumina.com

Illumina Customer Support Telephone Numbers

Region	Toll Free	Regional
North America	+1.800.809.4566	
Australia	+1.800.775.688	
Austria	+43 800006249	+43 19286540
Belgium	+32 80077160	+32 34002973
China	400.635.9898	
Denmark	+45 80820183	+45 89871156
Finland	+358 800918363	+358 974790110
France	+33 805102193	+33 170770446
Germany	+49 8001014940	+49 8938035677
Hong Kong	800960230	
Ireland	+353 1800936608	+353 016950506
Italy	+39 800985513	+39 236003759
Japan	0800.111.5011	
Netherlands	+31 8000222493	+31 207132960
New Zealand	0800.451.650	
Norway	+47 800 16836	+47 21939693
Singapore	+1.800.579.2745	
Spain	+34 911899417	+34 800300143
Sweden	+46 850619671	+46 200883979
Switzerland	+41 565800000	+41 800200442
Taiwan	00806651752	
United Kingdom	+44 8000126019	+44 2073057197
Other countries	+44.1799.534000	

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Document # 100000037191 v00

[This page intentionally left blank]

Illumina
5200 Illumina Way
San Diego, California 92122 U.S.A.
+1.800.809.ILMN (4566)
+1.858.202.4566 (outside North America)
techsupport@illumina.com

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
Reagent kit Praxis™ Extended RAS Panel for qualitative determination of the status of specific mutations genes of RAS (KRAS (exons 2, 3, and 4) and NRAS (exons 2, 3, and 4) by sequencing a new generation in a human genomic DNA sample isolated from colorectal cancer tissue samples fixed in formalin and enclosed in paraffin blocks

Instruction for Use

Manufacturer:


Illumina, Inc

Address:

5200 Illumina Way, San Diego, California 92122, USA


Approved

On behalf of Illumina, Inc.

Aurelia Brownridge
Sr. Regulatory Affairs Specialist
Illumina, Inc.

Date May 14, 2020

ACKNOWLEDGMENT

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California
County of San Diego)

On May 14 2020 before me, Mary Agostino, Notary Public
(insert name and title of the officer)

personally appeared Alexia Browne,
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are
subscribed to the within instrument and acknowledged to me that he/she/they executed the same in
his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the
person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing
paragraph is true and correct.

WITNESS my hand and official seal.

Signature Mary Agostino (Seal)

Praxis™ Extended RAS Panel

FOR *IN VITRO* DIAGNOSTIC USE

Catalog # 20012431: 2 Runs, up to 10 samples plus 2 controls per run/20 samples per kit

Intended Use

The Praxis™ Extended RAS Panel is a qualitative *in vitro* diagnostic test using targeted high throughput parallel sequencing for the detection of 56 specific mutations in RAS genes [*KRAS* (exons 2, 3, and 4) and *NRAS* (exons 2, 3, and 4)] in DNA extracted from formalin-fixed, paraffin-embedded (FFPE) colorectal cancer (CRC) tissue samples. The test is intended to be used on the Illumina MiSeqDx® instrument.

Summary and Explanation of the Assay


Clinical Description

Mutations in the RAS oncogenes are frequently found in human cancers. These mutations activate the mitogen-activated protein kinase pathway contributing to the malignant phenotype. Historically, *KRAS* exon 2 (codon 12 and 13) mutations were an established predictive biomarker of resistance to anti-epidermal growth factor receptor (EGFR) therapy for metastatic colorectal cancer (mCRC). Multiple clinical trials have since established the value of extended RAS testing to include exon 2 (codons 12 and 13), exon 3 (codons 59 and 61), and exon 4 (codons 117 and 146) of either *KRAS*, or *NRAS*.¹ Overall, extended RAS mutations (*KRAS* exons 3, 4 and *NRAS* exons 2, 3, 4) account for 11% of patients with mCRC.¹ Taken together, mutations in *NRAS* and *KRAS* may affect up to approximately 53% of patients.^{1,2,3}

RAS Mutations in Panel

The Praxis Extended RAS Panel is an *in vitro* diagnostic test used in the detection of specific mutations in exons 2, 3, and 4 of *KRAS* and *NRAS*. Figure 1 depicts the general locations of the mutations within the exons of the genes. Codon changes with red arrows show location of mutations targeted in the assay.

Figure 1 RAS Genes (*KRAS* and *NRAS*) and Targeted Sequencing

The assay is designed specifically for the detection of the mutations from the *KRAS* and *NRAS* genes listed in Table 1. The assay simultaneously determines the presence or absence of each of the 56 mutations in every sequenced sample.

Praxis Extended RAS Panel

Praxis Extended RAS Panel Summary of Mutations*

Chromosome	Gene	Exon	Amino Acid	Wild Type > Mutation
chr12	KRAS	2	p.Gly13Val	c.38_39GC>TT
chr12	KRAS	2	p.Gly13Glu	c.38_39GC>AA
chr12	KRAS	2	p.Gly13Asp	c.38_39GC>AT
chr12	KRAS	2	p.Gly13Asp	c.38C>A
chr12	KRAS	2	p.Gly13Cys	c.37C>T
chr12	KRAS	2	p.Gly13Arg	c.37C>C
chr12	KRAS	2	p.Gly12Trp	c.34_36GGT>TGG
chr12	KRAS	2	p.Gly12Ala	c.35C>C
chr12	KRAS	2	p.Gly12Asp	c.35C>A
chr12	KRAS	2	p.Gly12Phe	c.34_35GG>TT
chr12	KRAS	2	p.Gly12Asn	c.34_35GG>AA
chr12	KRAS	2	p.Gly12Val	c.35C>T
chr12	KRAS	2	p.Gly12Ser	c.34C>A
chr12	KRAS	2	p.Gly12Cys	c.34C>T
chr12	KRAS	2	p.Gly12Arg	c.34C>C
chr12	KRAS	3	p.Gln61His	c.183A>C
chr12	KRAS	3	p.Gln61His	c.183A>T
chr12	KRAS	3	p.Gln61Leu	c.182A>T
chr12	KRAS	3	p.Gln61Arg	c.182A>G
chr12	KRAS	3	p.Gln61Lys	c.181C>A
chr12	KRAS	3	p.Gln61Glu	c.181C>G
chr12	KRAS	3	p.Ala59Gly	c.176C>G
chr12	KRAS	3	p.Ala59Thr	c.175G>A
chr12	KRAS	4	p.Ala146Val	c.437C>T
chr12	KRAS	4	p.Ala146Thr	c.436G>A
chr12	KRAS	4	p.Ala146Pro	c.436G>C
chr12	KRAS	4	p.Lys117Asn	c.351A>T
chr12	KRAS	4	p.Lys117Asn	c.351A>C
chr1	NRAS	2	p.Gly13Glu	c.38_39GT>AG
chr1	NRAS	2	p.Gly13Glu	c.38_39GT>AA
chr1	NRAS	2	p.Gly13Asp	c.38C>A
chr1	NRAS	2	p.Gly13Val	c.38C>T
chr1	NRAS	2	p.Gly13Cys	c.37C>T
chr1	NRAS	2	p.Gly13Arg	c.37C>C
chr1	NRAS	2	p.Gly12Trp	c.34_36GGT>TGG
chr1	NRAS	2	p.Gly12Ala	c.35C>C
chr1	NRAS	2	p.Gly12Val	c.35C>T
chr1	NRAS	2	p.Gly12Asp	c.35C>A
chr1	NRAS	2	p.Gly12Phe	c.34_35GG>TT

Chromosome	Gene	Exon	Amino Acid	Wild Type > Mutation
chr1	NRAS	2	p.Gly12Asn	c.34_35GG>AA
chr1	NRAS	2	p.Gly12Ser	c.34G>A
chr1	NRAS	2	p.Gly12Cys	c.34G>T
chr1	NRAS	2	p.Gly12Arg	c.34G>C
chr1	NRAS	3	p.Gln61His	c.183A>T
chr1	NRAS	3	p.Gln61His	c.183A>C
chr1	NRAS	3	p.Gln61Arg	c.182A>G
chr1	NRAS	3	p.Gln61Leu	c.182A>T
chr1	NRAS	3	p.Gln61Lys	c.181C>A
chr1	NRAS	3	p.Gln61Glu	c.181C>G
chr1	NRAS	3	p.Ala59Gly	c.176C>G
chr1	NRAS	3	p.Ala59Thr	c.175G>A
chr1	NRAS	4	p.Ala146Val	c.437C>T
chr1	NRAS	4	p.Ala146Thr	c.436G>A
chr1	NRAS	4	p.Ala146Pro	c.436G>C
chr1	NRAS	4	p.Lys117Asn	c.351G>C
chr1	NRAS	4	p.Lys117Asn	c.351G>T

* KRAS proto-oncogene, GTPase [Homo sapiens (human)] (Gene ID: 3845) - NM_004985.4, NP_004976.2.

* NRAS neuroblastoma RAS viral oncogene homolog [Homo sapiens (human)] (Gene ID: 4893) - NM_002524.4, NP_002515.1.

Testing for KRAS and NRAS mutations is recommended by the NCCN Guidelines® for all patients with mCRC. Only patients who have wild type RAS are eligible for anti-EGFR therapy⁴.


Principles of Procedure

The Illumina Praxis Extended RAS Panel involves three main procedures. The first is to qualify the DNA sample that can be used for the test (Figure 2). The second is to manually prepare the samples for sequencing, which is called library preparation. Library preparation consists of four key steps: Hybridization, Extension-Ligation, PCR Amplification, and Library Normalization (Figure 3). The third procedure is to sequence the prepared sample using SBS (sequencing by synthesis) chemistry on the MiSeqDx.

Qualification of FFPE Samples


The Praxis Extended RAS Panel workflow starts with DNA qualification step which is performed to determine if the DNA is of sufficient quality for use with the Praxis Extended RAS Panel. Genomic DNA samples are qualified by measuring amplifiability relative to a control DNA template.

Figure 2 Amplification Graph

Library Preparation

Figure 3 Library Preparation Process

- A Hybridization**—A pool of upstream and downstream oligonucleotides specific to the targeted regions in the Praxis Extended RAS Panel is hybridized to the sample genomic DNA. Two oligonucleotide pools (A and B) independently test each strand of DNA. At the end of this process, a three-step wash procedure with a filter

capable of size selection removes unbound oligonucleotides from the genomic DNA.

- B **Extension-Ligation**—The second step, Extension-Ligation, connects the hybridized upstream and downstream oligonucleotides. A DNA polymerase extends from the upstream oligonucleotides through the targeted region, followed by ligation to the 5' end of the downstream oligonucleotide using a DNA ligase. The result is the formation of products that contain the RAS-specific oligonucleotides flanked by sequences required for amplification.
- C **PCR Amplification**—The third step, PCR Amplification, amplifies the extension-ligation products using primers that add index sequences for sample multiplexing, as well as common adapters required for cluster generation on the MiSeqDx. At the end of this process, a PCR clean-up procedure purifies the PCR products (referred to as a library).
- D **Library Normalization**—The final step, Library Normalization, normalizes the quantity of each library to ensure more equal library representation in the final pooled library. At the end of this process, the pooled library is loaded onto the MiSeqDx for sequencing using SBS chemistry.

Sequencing

SBS chemistry uses a reversible-terminator method to detect single nucleotide bases as they are incorporated into growing DNA strands. During each sequencing cycle, a single fluorescently-labeled deoxynucleotide triphosphate (dNTP) is added to the nucleic acid chain. The nucleotide label serves as a terminator for polymerization, so after each dNTP incorporation, the fluorescent dye is imaged to identify the base and then enzymatically cleaved to allow incorporation of the next nucleotide. Because all four- reversible terminator-bound dNTPs (A, G, T, C) are present as single, separate molecules, natural competition minimizes incorporation bias. Base calls are made directly from signal intensity measurements during each sequencing cycle. The result is base-by-base sequencing.

Data Analysis

The Praxis Extended RAS Panel module of the Local Run Manager software processes base calls generated during primary analysis and produces information about each sample, called secondary analysis. Secondary analysis includes demultiplexing, FASTQ file generation, alignment, variant calling, and report generation.

- ▶ **Demultiplexing**— This is the first step in secondary analysis. Demultiplexing separates data from pooled samples based on the unique sequence indexes that were added during the PCR amplification step.
- ▶ **FASTQ File Generation**—After demultiplexing, the Praxis Extended RAS Panel module generates intermediate files in the FASTQ format, which is a text format used to represent sequences. FASTQ files contain the reads for each sample and the quality scores, excluding reads from any clusters that did not pass filter.
- ▶ **Alignment**—Alignment compares sequences against the reference to identify a relationship between the sequences and assigns a score based on regions of similarity. Aligned reads are written to files in BAM format. Local Run Manager uses a banded Smith-Waterman algorithm that performs local sequence alignments to determine similar regions between two sequences.
- ▶ **Variant Calling**—Developed by Illumina, the somatic variant caller identifies mutations present at low frequency in the DNA sample. In order to minimize false positives due to FFPE artifacts, mutations are first called for each pool separately. Then, mutations from the two pools are compared and combined into a single output file. If found, only mutations listed in Table 1 are reported. For more information, see the *Praxis Extended RAS Panel Reference Guide (document # 1000000037191)*.
- ▶ **Report Generation**—Results are accessible through Local Run Manager in the Samples and Results tab for the run, and the following simple-text output file is generated:
 - ▶ **Extended RAS Panel Report**—Details whether each sample contains any of the mutations listed in the reference panel.

Limitations of the Procedure

- 1 For *in vitro* diagnostic use.
- 2 In this procedure, the **Extended RAS Panel Report** only determines the presence or absence of the NRAS and KRAS mutations listed in Table 1.

- 3 In a study assessing 685 banked specimens, 31% failed to meet the QC requirement of $dCq \leq 5.0$.
- 4 Neoplastic content is required to be $\geq 50\%$ in order for a FFPE tissue to be processed for sequencing.

Reagents

Reagents Provided

Illumina provides the following reagents for the Illumina Praxis Extended RAS Panel (Table 2–Table 10). Catalog # 20012431 has been configured for 2 runs with a maximum of 10 samples plus 2 controls per run.

Praxis Extended RAS Panel, Box 1

Table 2 Box 1A Pre-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Hybridization Buffer	1 tube	4.32 ml	Buffered aqueous solution containing salts and formamide	-25°C to -15°C
Extension-Ligation Mix	1 tube	4.8 ml	Buffered aqueous solution containing proprietary blend of DNA polymerases, DNA ligase, and dNTPs	-25°C to -15°C
Index Primers A-H (A501-A508)	1 tube per primer	192 μ l	PCR primers with index sequences and sequencing adapters	-25°C to -15°C
Index Primers 1-12 (A701-A712)	1 tube per primer	128 μ l	PCR primers with index sequences and sequencing adapters	-25°C to -15°C
PCR Polymerase	1 tube	56 μ l	Proprietary DNA polymerase	-25°C to -15°C
PCR Master Mix	1 tube	2.8 ml	Buffered aqueous solution containing salts and dNTPs	-25°C to -15°C

Table 3 Box 1B Post-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Library Normalization Diluent	1 tube	4.6 ml	Buffered aqueous solution containing salts, 2-Mercaptoethanol, and formamide	-25°C to -15°C
Library Dilution Buffer	2 tubes	4.5 ml	Buffered aqueous solution	-25°C to -15°C
PhiX Internal Control	1 tube	10 μ l	Buffered aqueous solution containing PhiX genomic DNA	-25°C to -15°C

Praxis Extended RAS Panel, Box 2

Table 4 Box 2 Post-Amp Reagents

Component	Quantity	Contents	Storage
MiSeq™ Dx SBS Cartridge v3	2 cartridges	Single-use cartridge that contains cluster generation and sequencing reagents for use with the MiSeqDx, including formamide, 2-Mercaptoethanol, and < 2% DMSO	-25°C to -15°C

Praxis Extended RAS Panel, Box 3

Table 3 Box 3A Pre-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Stringent Wash Buffer	1 bottle	24 ml	Buffered aqueous solution containing salts, 2-Mercaptoethanol and formamide	2°C to 8°C
Universal Wash Buffer	1 tube	4.8 ml	Buffered aqueous solution containing salts	2°C to 8°C

Table 4 Box 3B Post-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
PCR Clean-Up Beads	1 tube	5 ml	Buffered aqueous solution containing solid phase paramagnetic beads and polyethylene glycol	2°C to 8°C
Library Normalization Wash	2 tubes	4.8 ml	Buffered aqueous solution containing salts, 2-Mercaptoethanol and formamide	2°C to 8°C
Library Beads	1 tube	1.2 ml	Buffered aqueous solution containing solid phase paramagnetic beads	2°C to 8°C
MiSeq™ Dx Flow Cell	2 containers	1 flow cell	Glass substrate with covalently bound oligonucleotides	2°C to 8°C

Praxis Extended RAS Panel, Box 4

Table 5 Box 4 Post-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
MiSeq™ Dx SBS Solution (PR2)	2 bottles	500 ml	Buffered aqueous solution	2°C to 8°C

Praxis Extended RAS Panel, Box 5

Table 6 Box 5 Pre-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
TruSeq™ Custom Amplicon Dx Filter Plate	2 plates	N/A	Polypropylene microtiter plate with a modified polyethersulfone membrane	15°C to 30°C

Table 7 Box 5 Post-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Elution Buffer	1 tube	4.8 ml	Buffered aqueous solution	15°C to 30°C
Library Storage Buffer	1 tube	3.5 ml	Buffered aqueous solution	15°C to 30°C

Praxis Extended RAS Panel, Box 6

Table 1 Box 6 Pre-Amp Reagents

Component	Quantity	Fill Volume	Active Ingredients	Storage
Oligo Pool A—Extended RAS	1 tube	150 µl	Buffered aqueous solution containing oligonucleotides targeting the RAS genes	-25°C to -15°C
Oligo Pool B—Extended RAS	1 tube	150 µl	Buffered aqueous solution containing oligonucleotides targeting the RAS genes	-25°C to -15°C
Quality Control (QC) Primers	4 tubes	75 µl	Buffered aqueous solution containing oligonucleotides for DNA sample qualification	-25°C to -15°C
Quality Control (QC) Template	1 tube	220 µl	Buffered aqueous solution containing genomic DNA for sample qualification	-25°C to -15°C
qPCR Master Mix	2 tubes	1 ml	Buffered aqueous solution containing qPCR reaction components	-25°C to -15°C
Panel Control—Extended RAS	1 tube	8 µl	Buffered aqueous solution containing genomic DNA for assay positive control	-25°C to -15°C

Reagents Required, Not Provided

Pre-Amp Reagents

- ✦ 10 N NaOH (prepare from tablets or use a standard solution)
- ✦ TE Buffer
- ✦ RNase/DNase-free water

Post-Amp Reagents

- ✦ 10 N NaOH (prepare from tablets or use a standard solution)
- ✦ Ethanol, 200 proof for molecular biology
- ✦ TE Buffer
- ✦ RNase/DNase-free water

Storage and Handling

- 1 Room temperature is defined as 15°C to 30°C.
- 2 Reagents are provided in sufficient quantity for 2 separate uses of a maximum of 12 samples/24 libraries each, including controls. Between uses, reagents should be recapped and returned to their original storage temperatures, unless otherwise indicated as single use only.
- 3 The following reagents are shipped frozen and are stable when stored at -25°C to -15°C until the specified expiration date.
 - ✦ Oligo Pool A—Extended RAS
 - ✦ Oligo Pool B—Extended RAS
 - ✦ Hybridization Buffer
 - ✦ Extension-Ligation Mix
 - ✦ Index Primers A-H (A501-A508)
 - ✦ Index Primers 1-12 (A701-A712)
 - ✦ PCR Polymerase
 - ✦ PCR Master Mix
 - ✦ Library Normalization Diluent
 - ✦ Library Dilution Buffer
 - ✦ PhiX Internal Control
 - ✦ MiSeq™ Dx SBS Cartridge v3
 - ✦ QC Primers
 - ✦ QC Template

- ▶ Panel Control—Extended RAS
- ▶ qPCR Master Mix

Except for the reagent cartridge, the reagents above are stable for a maximum of 5 freeze/thaw cycles that occur before the specified expiration date.

Index primers and the MiSeq™ Dx SBS Cartridge v3 are single-use only.

CAUTION

Do not store reagents in a frost-free storage unit or store in door compartments.

Do not refreeze the reagent cartridge after it has been thawed. It can be stored up to 6 hours at 2°C to 8°C.

- 4 The following reagents are shipped refrigerated and are stable when stored at 2°C to 8°C until the specified expiration date.
 - ▶ Stringent Wash Buffer
 - ▶ Universal Wash Buffer
 - ▶ PCR Clean-Up Beads
 - ▶ Library Beads
 - ▶ Library Normalization Wash
 - ▶ MiSeq™ Dx SBS Solution (PR2)
 - ▶ MiSeq™ Dx Flow Cell

SBS Solution and Flow Cell are single-use only.
- 5 The following reagents are shipped ambient and are stable when stored at room temperature until the specified expiration date:
 - ▶ Elution Buffer
 - ▶ TruSeq™ Custom Amlicon Dx Filter Plate
 - ▶ Library Storage Buffer

Filter plate is single-use only.
- 6 Changes in the physical appearance of the reagents provided can indicate deterioration of the materials. If changes in the physical appearance occur (e.g., obvious changes in reagent color or cloudiness apparent with microbial contamination), do not use the reagents.
- 7 The Hybridization Buffer, Stringent Wash Buffer, and Library Normalization Diluent reagents might form visible precipitates or crystals. Before use, vortex vigorously, and then visually inspect to ensure that no precipitates are present.
- 8 Adhere to the following best practices when handling PCR Clean-Up Beads and Library Beads:
 - ▶ The beads should never be frozen.
 - ▶ Allow the beads to reach room temperature.
 - ▶ Immediately before use, vortex the beads until well-suspended and color appears homogeneous.
 - ▶ Thoroughly mix the sample after the beads are added by pipetting up and down 10 times. Alternatively, use a shaker to thoroughly mix samples.
 - ▶ Incubate the bead/sample mixture at room temperature for the entire duration indicated.
 - ▶ Follow instructions when using the magnetic stand. Wait for the solution to clear before aspirating. Keep the plate on the magnetic stand when slowly aspirating the supernatant, taking care not to disturb the separated beads.
- 9 The PCR amplification plate can remain on the thermal cycler overnight, or it can be stored at 2°C to 8°C for up to 2 days. Before storing the plate at 2°C to 8°C, seal the plate well.

CAUTION

Do not freeze the Library Beads or mix with the Library Normalization Diluent reagent if not used immediately.

- 10 The diluted amplicon library can be stored at -25°C to -15°C for up to 30 days.

Equipment and Materials

Equipment and Materials Required, Sold Separately


- 1 MiSeq Dx Instrument, Catalog # DX-410-1001

- 2 **TruSeq® Index Plate Fixture Kit**, Catalog # FC-130-1005
- 3 **TruSeq Index Plate Fixture & Collar Kit**, Catalog # FC-130-1007

Equipment and Materials Required, Not Provided

Pre-Amp Equipment and Materials


- 1 **DNA Extraction and Purification Reagents**—See *DNA Extraction* on page 12.
- 2 **Heat Block**—One heat block for a 96-well PCR plate is required. The heat block must meet the following performance specifications. Heat blocks with heated lids are required for use.
 - ▶ Temperature Range: Ambient +5°C to 99°C
 - ▶ Temperature Regulation: $\pm 0.1^\circ\text{C}$ at 37°C ; $\pm 0.4^\circ\text{C}$ at 60°C
- 3 **Sample Incubator**—One incubator (hybridization oven) is required. The incubator must meet the following performance specifications.
 - ▶ Temperature Range: 10°C to 100°C
 - ▶ Temperature Regulation: $\pm 0.2^\circ\text{C}$
- 4 **Tabletop Centrifuge**—One tabletop centrifuge is required. (A separate centrifuge is required in the post-amp area.) Any plate centrifuge that attains the designated speeds of the protocol (280 to $2400 \times g$) is acceptable.
- 5 **Microcentrifuge**—One microcentrifuge for microcentrifuge tubes.
- 6 **Precision Pipettes**—One set of precision pipettes is required. (A separate set is required in the post-amp area.) The use of precision pipettes is required to ensure accurate reagent and sample delivery. If they are calibrated regularly and are accurate within 5% of stated volume, single-channel or multi-channel pipettes can be used.
- 7 **Sealing Rollers or Sealing Wedge**—A sealing roller or sealing wedge are required for sealing microplates.
- 8 **Consumables**—The following consumables are required.
 - ▶ 96-well or 384-well optical PCR plates and seals
 - ▶ 96-well skirted PCR plates, 0.2 ml, polypropylene, or equivalent
 - ▶ 96-well storage plates, 0.8 ml (MIDI plates)

NOTE

Make sure that the 96-well plate is fit compatible with the magnetic stand.

- ▶ PCR eight-tube strips
- ▶ Solution basin, PVC, RNase/DNase-free (trough)
- ▶ Thermal microplate sealer

NOTE

If a thermal microplate sealer is not available, adhesive aluminum foil seal can be used.

- ▶ Adhesive aluminum foil seal
- ▶ Appropriate PCR plate seals
- ▶ Aerosol resistant pipette tips

Post-Amp Equipment and Materials

- 1 **96-Well or 384-Well Real-Time PCR Detection System**—One Real-Time PCR Detection system is required. The system must meet the following performance specifications. In addition, the system must have a heated lid and support a 96-well or 384-well format.
 - ▶ Minimum Temperature Range: 37°C to 95°C
 - ▶ Accuracy/Uniformity: $\pm 0.5^\circ\text{C}$
 - ▶ Dye: FAM
- 2 **Thermal Cycler**—One thermal cycler is required. The thermal cycler must have a heated lid and meet the following performance specifications:
 - ▶ Minimum Temperature Control Range: 4°C to 99°C
 - ▶ Control Accuracy: $\pm 0.3^\circ\text{C}$
- 3 **Optical Compression Pad**—An optical compression pad is recommended to provide a proper thermal seal between the PCR thermal cycler and reaction plates when using adhesive plate seals.

- 4 **Microplate Shaker**—One microplate shaker is required in the post-amp lab area. The plate shaker must meet the following performance specifications:
 - ▶ Maximum Mixing Speed: 3000 rpm
 - ▶ Mixing Speed Range: 200 to 3000 rpm
 - 5 **Tabletop Centrifuge**—One tabletop centrifuge is required. (A separate centrifuge is required in the pre-amp area.) Any plate centrifuge that attains the designated speeds of the protocol (280 to 2400 × g) is acceptable.
 - 6 **Heat Block**—One heat block for tubes is required. The heat block must meet the following performance specifications.
 - ▶ Temperature Range: Ambient +5°C to 99°C
 - ▶ Temperature Regulation: ±0.1°C at 37°C; ±0.4°C at 60°C
 - 7 **Magnetic Stand**—One magnetic stand for a 96-well plate is required. Better performance is seen when the magnets are on the side of the stand and not on the bottom.
 - 8 **Precision Pipettes**—One set of precision pipettes is required. (A separate set is required in the pre-amp area.) The use of precision pipettes is required to ensure accurate reagent and sample delivery. If they are calibrated regularly and are accurate within 5% of stated volume, single-channel or multi-channel pipettes can be used.
 - 9 **Gel Electrophoresis Supplies**—Gel electrophoresis supplies and apparatus are required along with an appropriate staining method to visualize PCR products in the gel.
 - 10 **Consumables**—The following consumables are required.
 - ▶ 96-well skirted PCR plates, 0.2 ml, polypropylene, or equivalent
 - ▶ 96-well storage plates, 0.8 ml (MIDI plates)
-  **NOTE**
Make sure that the 96-well plate is fit compatible with the magnetic stand.
- ▶ 2-4% TBE Agarose gel
 - ▶ 100 bp DNA molecular weight marker
 - ▶ DNA loading dye
 - ▶ Conical tubes, 15 ml
 - ▶ Microcentrifuge tubes
 - ▶ PCR eight-tube strips
 - ▶ Solution basins, PVC, RNase/DNase-free (trough)
 - ▶ Adhesive plate seals
 - ▶ Aerosol resistant pipette tips

Specimen Collection, Transport, and Storage


Tissue Assessment

The Praxis Extended RAS Panel utilizes DNA samples extracted from formalin-fixed paraffin-embedded (FFPE) colorectal cancer (CRC) tissue samples.

To prepare FFPE CRC tissue samples for DNA extraction:

- 1 Scrape excess paraffin away from the tissue using a fresh, sterile scalpel. Use serial sections for DNA extraction.
- 2 Use sections with ≥ 50% neoplastic content by area for processing without macrodissection.
- 3 For sections that are less than 50% neoplastic content by area, enrich neoplastic content by macrodissecting the sections to obtain ≥ 50% neoplastic content by area. In order to maximize total tissue area captured per slide, macrodissecting to obtain 50% neoplastic content by area is recommended.

Figure 4 Maximizing Total Tissue Area Captured per Slide During Macrodissection

20% neoplastic content by area.

Macrodissect to obtain $\geq 50\%$ neoplastic content by area, while maximizing total tissue area captured per slide.

- 4 Cumulative tissue area is the total tissue area per section times the total number of $5\ \mu\text{m}$ sections. For optimal results, recommended cumulative tissue area using $5\ \mu\text{m}$ sections is $240\ \text{mm}^2$ or greater with $\geq 50\%$ neoplastic content by area. For optimal yield of amplifiable DNA, use at least $8 \times 5\ \mu\text{m}$ sections or equivalent alternate thickness serial sections.

**NOTE**


If macrodissection was performed, the macrodissected tissue area should be used in the determination of cumulative tissue area.

Recommended Cumulative Tissue Area is $\geq 240\ \text{mm}^2$

**Per Slide as Required:**

Macrodissect to obtain $\geq 50\%$ neoplastic content by area, while maximizing total tissue area captured per slide.

**Recommended Cumulative Tissue Area:**

For $5\ \mu\text{m}$ sections, $\geq 240\ \text{mm}^2$ cumulative tissue area with $\geq 50\%$ neoplastic content by area. For optimal DNA yield, use at least eight slides.


- 5 Minimum cumulative tissue area using $5\ \mu\text{m}$ sections is $80\ \text{mm}^2$ with $\geq 50\%$ neoplastic content by area. Equivalent cumulative tissue area for alternate thickness serial sections can also be used.

Minimum Cumulative Tissue Area is $80\ \text{mm}^2$

**Per Slide as Required:**

Macrodissect to obtain $\geq 50\%$ neoplastic content by area, while maximizing total tissue area captured per slide.

**Minimum Cumulative Tissue Area:**

For $5\ \mu\text{m}$ sections, $80\ \text{mm}^2$ cumulative tissue area with $\geq 50\%$ neoplastic content by area.

- 6 Extracted DNA samples giving a delta Cq value ≤ 5.00 can be used for the Praxis Extended RAS Panel.
7 [Optional] If extracted DNA samples do not meet the delta Cq threshold, additional sections can be used for extraction if available.

DNA Extraction

- ▶ Illumina recommends column-based DNA extraction kits, using double the amount of Proteinase K, overnight Proteinase K incubations with agitation, and final elutions in at least $30\ \mu\text{l}$ volume.
- ▶ Bead-based extraction methods and methods using only lysis of crude cell extracts are not recommended for use with this assay.

- ▶ Illumina recommends storing extracted genomic DNA between 2°C and 8°C for a maximum of 28 days or storing frozen between -15°C to -25°C for a maximum of 161 days. Genomic DNA stored frozen is stable for a maximum of 2 freeze/thaw cycles.

Warnings and Precautions

CAUTION

Federal law restricts this device to sale by or on the order of a physician or other practitioner licensed by the law of the State in which he/she practices, to use or order the use of the device.

- 1 Some components of this assay contain potentially hazardous chemicals. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Wear protective equipment, including eye protection, gloves, and laboratory coat appropriate for risk of exposure. Handle used reagents as chemical waste and discard in accordance with applicable regional, national, and local laws and regulations. For additional environmental, health, and safety information, see the SDS at support.illumina.com/sds.html. (See *Reagents* on page 6 for more information.)
- 2 Handle all specimens as if they are potentially infectious agents.
- 3 Failure to follow the procedures as outlined may result in erroneous results or significant reduction in sample quality.
- 4 Use routine laboratory precautions. Do not pipette by mouth. Do not eat, drink, or smoke in designated work areas. Wear disposable gloves and laboratory coats when handling specimens and assay reagents. Wash hands thoroughly after handling specimens and assay reagents.
- 5 Do not use any assay components beyond their stated expiration date on the assay carton label. Do not interchange assay components from different assay lots. Assay lots are identified on the assay carton label.
- 6 Store the assay components at the specified temperature in designated pre-amplification and post-amplification areas.
- 7 Repeated freeze-thaw cycles (up to 5) of the Box 1 and Box 6 components will not compromise the integrity of the assay.
- 8 Index primers are single use only and should be discarded after use.
- 9 To prevent sample or reagent degradation, ensure that all sodium hypochlorite vapors have fully dissipated before starting the protocol.
- 10 Proper laboratory practices and good laboratory hygiene is required to prevent PCR products from contaminating reagents, instrumentation, and genomic DNA samples. PCR contamination may cause inaccurate and unreliable results.
- 11 To prevent contamination, ensure that pre-amplification and post-amplification areas have dedicated equipment (e.g., pipettes, pipette tips, vortexer, and centrifuge).
- 12 Avoid cross contamination. Use fresh pipette tips between samples and between dispensing reagents. Mix samples with a pipette and centrifuge the plate when indicated. Do not vortex the plates. Using aerosol-resistant tips reduces the risk of amplicon carry-over and sample-to-sample cross-contamination.
- 13 Index-sample pairing must match the printed plate layout exactly. Mismatches between the printed plate layout and the sample will result in loss of positive sample identification and incorrect result reporting.
- 14 Always prepare fresh 80% ethanol for wash steps. Ethanol can absorb water from the air, impacting results.
- 15 Ensure that all ethanol is removed from the bottom of the wells during wash steps. Residual ethanol may impact results.
- 16 Adhere to the specified drying time following the magnetic stand step to ensure complete evaporation. Residual ethanol can impact the performance of subsequent reactions.
- 17 Do not mix the Oligo Pool A and Oligo Pool B and Hybridization Buffer for storage. When combined, the Oligo Pools become unstable, even when stored frozen.
- 18 The use of thermal cyclers with active cooling (e.g., Peltier, thermoelectric cooled) is not recommended for the hybridization step. The passive cooling step is critical for proper hybridization.
- 19 Always add PCR Polymerase to the PCR Master Mix immediately before use. Never store the combined working solution.
- 20 During the library normalization step, it is critical to resuspend the library bead pellet. This is essential for achieving consistent cluster density on the MiSeq™ Dx flow cell.

- 21 Adhere to the specified incubation times in the library normalization step. Improper incubation can affect library representation and cluster density.
- 22 Due to the number of plate transfers and subsequent potential for contamination, extreme care should be taken to ensure that well contents remain fully in the well. Do not splash contents.
- 23 DNA samples giving a delta C_q value less than or equal to 5.00 can be used for the Praxis Extended RAS Panel. Illumina does not support the use of gDNA samples giving a delta C_q value of greater than 5.00 with the Praxis Extended RAS Panel.

Index Primer Sequences

Sequences for Index Primers A-H (A501-A508)

Index Primer	Sequence
Index Primer A (A501)	TGAACCTT
Index Primer B (A502)	TGCTAAGT
Index Primer C (A503)	TGTTCTCT
Index Primer D (A504)	TAAGACAC
Index Primer E (A505)	CTAATCGA
Index Primer F (A506)	CTAGAACA
Index Primer G (A507)	TAAGTTCC
Index Primer H (A508)	TAGACCTA

Sequences for Index Primers 1-12 (A701-A712)

Index Primer	Sequence
Index Primer 1 (A701)	ATCACGAC
Index Primer 2 (A702)	ACAGTGGT
Index Primer 3 (A703)	CAGATCCA
Index Primer 4 (A704)	ACAAACGG
Index Primer 5 (A705)	ACCCAGCA
Index Primer 6 (A706)	AACCCCTC
Index Primer 7 (A707)	CCCAACCT
Index Primer 8 (A708)	CACCACAC
Index Primer 9 (A709)	GAAACCCA
Index Primer 10 (A710)	TGTGACCA
Index Primer 11 (A711)	AGGGTCAA
Index Primer 12 (A712)	AGGAGTGG

Instructions for Use

Qualification of DNA Extracted from FFPE Samples

Preparation

- 1 Bring the QC Template, QC Primers, qPCR Master Mix, and genomic DNA to room temperature.
- 2 Place thawed tubes on ice.

NOTE

If using QC Template for the first time, aliquot 7 µl of QC Template into different PCR strip tubes for long-term storage to avoid freeze-thawing.

Procedure

- 1 Determine the plate layout of the qPCR reaction.

**NOTE**

Samples, QC Template, and No Template Control (NTC) all must be tested in triplicate. See the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191) for more information on plate layouts.

- 2 Add 5 µl of QC Template to 495 µl of RNase/DNase-free water in a microcentrifuge tube.
- 3 Vortex the dilution to mix the sample.
- 4 Add 1.5 µl of extracted genomic DNA to 148.5 µl of RNase/DNase-free water in a PCR eight-tube strip to make a 100-fold dilution.
- 5 Using a P200, pipette up and down 10 times to mix the dilutions.
- 6 Add 150 µl of the diluted QC Template to an unused well in the PCR eight-tube strip consistent with the plate layout.
- 7 To another unused well, add 150 µl of nuclease-free water for use as NTC sample.
- 8 Prepare the qPCR Master Mix reaction as follows for a 384-well or 96-well format. The master mix contains extra volume.
 - ▶ 384-Well Format
 - ▶ qPCR Master Mix: 5.0 µl/well (205 µl for 36 reactions)
 - ▶ QC Primers: 0.8 µl/well (33 µl for 36 reactions)
 - ▶ RNase/DNase-free Water: 2.2 µl/well (90 µl for 36 reactions)
 - ▶ 96-Well Format
 - ▶ qPCR Master Mix: 10.0 µl/well (410 µl for 36 reactions)
 - ▶ QC Primers: 1.6 µl/well (66 µl for 36 reactions)
 - ▶ RNase/DNase-free Water: 4.4 µl/well (180 µl for 36 reactions)

**NOTE**

See the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191) for more information on preparing the qPCR Master Mix reaction.

- 9 Mix reaction mix gently but thoroughly.
- 10 Place the reaction mix on ice and protect it from light until use.
- 11 Aliquot master mix into a trough or eight-tube strip to aid in dispensing with multi-channel pipettor.
- 12 Add 8 µl (384-well format) or 16 µl (96-well format) of the master mix to each well of the qPCR plate. Take care to pipette accurately into the wells as small variations will affect the assay.
- 13 Add 2 µl (384-well format) or 4 µl (96-well format) of the QC Template dilution, the sample dilutions, or nuclease-free water to each well of the plate. Take care to pipette accurately into the wells as small variations will affect the assay.
- 14 Seal the plate with the appropriate seal for your instrument, taking care to avoid cross-contamination and to avoid smudging the surface of the lids. Verify that the seal is oriented correctly and uniformly sealed.
- 15 Place the plate on an adapter (if needed) and centrifuge the plate to 1000 × g for 1 minute.
- 16 Ensure that the seal is free of any liquid or dust, place the plate on the qPCR machine in the correct orientation, then close the lid and run the following qPCR thermal profile:
 - ▶ Hot Start
 - ▶ 50°C for 2 minutes
 - ▶ 95°C for 10 minutes
 - ▶ x40
 - ▶ 95°C for 30 seconds
 - ▶ 57°C for 30 seconds
 - ▶ 72°C for 30 seconds
 - ▶ Confirm that the instrument captures images after the 72°C step.

**NOTE**

The C_q threshold should be set to a value appropriate for your instrument that avoids inaccurate measurements due to background.

- 17 After the final step, the thermal cycler analyzes the samples. Ensure that the QC Template produces expected amplification curves. On most systems, the QC Template should amplify at a C_q of approximately 15–22 cycles. Exclude replicates from a triplicate group that are > 0.5 C_q different from the rest of the group. At least two of the three replicates less than or equal to 0.5 C_q apart must be included

in the final calculation. If these criteria are not met, the qualification process should be repeated for those samples.

**NOTE**

Four or more samples per 10-sample run having replicates removed indicate technical errors and the qualification process for all samples should be repeated.

**NOTE**

If little to no measurable amplification occurs for two or more replicates, the qualification process should be repeated for those samples.

- 18 Replicates exhibiting abnormal amplification curves should be excluded.

**NOTE**

See the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191) for more information on appropriate amplification curves.

- 19 Subtract the average Cq for the QC Template from the average Cq for each sample (Sample Average Cq – QC Template Average Cq) to yield the delta Cq values for each sample. Record the delta Cq values and the sample dilution factors. The sample dilution factors are based on the recommendations in the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191).
- 20 Make sure that amplification of the NTC occurs at least 10 cycles after QC Template amplification.

**NOTE**

Repeat entire qualification process for all samples if NTC fails.

SAFE STOPPING POINT

Illumina recommends storing extracted genomic DNA between 2°C to 8°C for a maximum of 28 days or storing frozen between –15°C to –25°C for a maximum of 161 days. Genomic DNA stored frozen is stable for a maximum of 2 freeze/thaw cycles.

Set Parameters

For more information on the Local Run Manager software, see the *Local Run Manager Software Reference Guide for MiSeqDx* (document # 1000000011880) and the *Local Run Manager Praxis Extended RAS Panel Analysis Module Workflow Guide* (document # 1000000011875).

- 1 Log in to Local Run Manager.
- 2 Click **Create Run**, and select **Extended RAS**.
- 3 Make sure that "Create Run Praxis Extended RAS" is visible in the top left corner of the screen.
- 4 Enter a run name that identifies the run from sequencing through analysis.
Use alphanumeric characters, spaces, underscores, or dashes.
- 5 [Optional] Enter a run description to help identify the run.
Use alphanumeric characters.

Specify Samples for the Run

Specify samples for the run using the following options:

- ▶ **Enter samples manually**—Use the table on the Create Run screen.
- ▶ **Import samples**—Navigate to an external file in a comma-separated values (*.csv) format.

Enter Samples Manually

- 1 Click the plate region A1 or E7 to enter sample information in the samples table.
Region A1 is equivalent to the upper left quadrant of wells in the Hyb plate layout, and includes samples for Oligo Pools A and B.
Region E7 is equivalent to the lower right quadrant of wells in the Hyb plate layout, and includes samples for Oligo Pools A and B.
- 2 In the Pool A table, enter a unique sample name in the Sample Name field for each well you use.
Use alphanumeric characters, dashes, or underscores.
The sample name automatically populates the corresponding well in the Pool B table.
Dual index adapters for each Index Read, specific to the well, display as you enter a sample name.

- 3 [Optional] Although the default layout is recommended, to change the position of control samples, use the Positive Control or Negative Control drop down menu.
- 4 Click the print  icon to display the plate layout.
- 5 Select **Print** to print the plate layout as a reference for preparing libraries.
- 6 Click **Close**.
- 7 Click **Save Run**.

Import Samples

- 1 Click the plate region A1 or E7.
- 2 Click the **Export**  icon for an example template.
- 3 Enter a unique sample name in the template for up to 10 samples and save the file.
Use alphanumeric characters, dashes, or underscores.

[Data]	Sample Names as a 4 * 3 matrix	
Sample1	Sample5	Sample9
Sample2	Sample6	Sample10
Sample3	Sample7	Positive_Control
Sample4	Sample8	Negative_Control

- 4 Click **Import Samples**.
- 5 Browse to the location of the sample information file in *.csv file format.
- 6 When finished, click the print  icon to display the plate layout.
- 7 Select **Print** to print the plate layout as a reference for preparing libraries.
- 8 Click **Close**.
- 9 Click **Save Run**.

Hybridization of Oligonucleotide Pool

Preparation

- 1 Determine genomic DNA samples to be used in the assay. Samples giving a delta Cq value less than or equal to 5.00 can be used for the Praxis Extended RAS Panel.
- 2 [Optional] If extracted DNA samples do not meet the delta Cq threshold, additional tissue sections can be used for DNA extraction if available.
- 3 Bring the Oligo Pool A, Oligo Pool B, Hybridization Buffer, genomic DNA samples, and Panel Control sample to room temperature.
- 4 Vortex the Oligo Pool A, Oligo Pool B, and Hybridization Buffer vigorously to make sure that all precipitates have dissolved, then briefly centrifuge the tubes to collect liquid.
- 5 Set a 96-well heat block to 95°C.
- 6 Pre-heat an incubator to 37°C.
- 7 Create the run and enter sample plate information in the Praxis Extended RAS Panel analysis module. Verify the location of positive and negative controls and then print sample plate layout.

Procedure

- 1 Set out a new 96-well PCR plate (hereafter referred to as **HYB** plate).

NOTE

See the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191) for more information on suggested **HYB** plate layout.

- 2 Add 10 µl of TE Buffer to each of the Negative Control (No Template Control) wells for Oligo Pool A and Oligo Pool B.
- 3 For each 12 sample/ 24 library preparation, add 2.5 µl of Panel Control to 22.5 µl of TE Buffer in a microcentrifuge tube. Mix the diluted Panel Control thoroughly.
- 4 Add 10 µl of diluted Panel Control prepared in the previous step to each of the two Positive Control wells for Oligo Pool A and Oligo Pool B.
- 5 Dilute genomic DNA samples (if necessary) according to the recommendations in Table 13.

Table 13 Delta C_q Dilution Table

Delta C _q	< -1.51	-1.50 to -0.51	-0.50 to 0.50	0.51 to 1.50	1.51 to 5.00	> 5.00
Dilution	4x	2x	No dilution	No dilution	No dilution	Not supported

- 6 Add 10 µl of each genomic DNA sample (diluted as necessary) to the appropriate wells for Oligo Pool A and Oligo Pool B.

**NOTE**

For information on diluting your samples and the suggested plate layout for correct Oligo Pool A and Oligo Pool B well selection, see the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191).

- 7 Add 5 µl of the Oligo Pool A to all Oligo Pool A wells containing genomic DNA, Panel Control, and Negative Control. Then, add 5 µl of the Oligo Pool B to all Oligo Pool B wells containing genomic DNA, Panel Control, and Negative Control.
- 8 Add 40 µl of Hybridization Buffer to each sample in the **HYB** plate. Gently pipette up and down 3–5 times to mix.

**NOTE**

Ensure any crystals or precipitate in Hybridization Buffer have dissolved.

- 9 Seal the **HYB** plate and centrifuge 1000 × g at room temperature for 1 minute.
- 10 Place the **HYB** plate in the pre-heated block at 95°C and incubate for 1 minute.
- 11 Change the temperature of the same heat block to 40°C and incubate for 3–18 hours.

**NOTE**

Moving the plate from the 95°C heat block to another pre-heated block set to 40°C will adversely affect hybridization.

Removal of Unbound Oligonucleotides

Preparation

- Bring Extension-Ligation Mix, Stringent Wash Buffer, and Universal Wash Buffer to room temperature, and then vortex briefly.
- Assemble the filter plate assembly (hereafter referred to as the **FPU**) in order from bottom to top: MIDI plate, adapter collar, filter plate, and lid.
- Pre-wash the filter plate membrane as follows:
 - Add 50 µl of Stringent Wash Buffer to each well.
 - Cover the filter plate with the lid and centrifuge at 2400 × g at room temperature for 5 minutes.
 - If the wash buffer does not drain completely, centrifuge the filter plate again at 2400 × g at room temperature for 5 minutes.

Procedure

- Remove the **HYB** plate from the heat block and centrifuge at 1000 × g at room temperature for 1 minute.
- Transfer the entire volume (approximately 65 µl) of each sample to the corresponding wells of the filter plate.
- Cover the filter plate with the lid and centrifuge at 2400 × g at room temperature for 5 minutes.

**NOTE**

Make sure that liquid has drained after centrifugation. Repeat centrifugation if necessary.

- Wash the filter plate as follows:
 - Add 50 µl of Stringent Wash Buffer to each sample well.
 - Cover the filter plate with the lid and centrifuge at 2400 × g at room temperature for 5 minutes.
 - If the wash buffer does not drain completely, centrifuge the filter plate again at 2400 × g at room temperature for 5 minutes.

- 5 Repeat the wash as described in the previous step.

**NOTE**

If the wash buffer does not drain completely, centrifuge again at $2400 \times g$ at room temperature until all liquid has gone through (an additional 5–10 minutes).

- 6 Discard all the flow-through (containing formamide), then reassemble the **FPU**.
- 7 Add 45 μ l of Universal Wash Buffer to each sample well.
- 8 Cover the filter plate with the lid and centrifuge at $2400 \times g$ at room temperature for 5 minutes.

**NOTE**

Make sure that all liquid has drained after centrifugation. Repeat centrifugation if necessary.

Extension-Ligation of Bound Oligonucleotides

Procedure

- 1 Add 45 μ l of Extension-Ligation Mix to each sample well of the filter plate.
- 2 Seal the filter plate with adhesive aluminum foil, and then cover with the lid.
- 3 Incubate the **FPU** in the pre-heated 37°C incubator for 45 minutes.
- 4 While the **FPU** plate is incubating, prepare the AMP (Amplification Plate) as described in the following section.

PCR Amplification

Preparation

- 1 Prepare fresh 0.05 N NaOH by adding 25 μ l of 10 N NaOH to 4975 μ l of sterile water.
- 2 Determine the index primers to be used according to the plate graphic printout from Praxis Extended RAS Panel analysis module.
- 3 Bring PCR Master Mix and the appropriate index primers to room temperature. Vortex each thawed tube to mix, and then briefly centrifuge the tubes.
- 4 Arrange the primers in a rack using the following arrangements:
 - a Arrange Index Primers A-H (A501-A508) primer tubes (white caps, clear solution) vertically, aligned with rows A through H.
 - b Arrange Index Primers 1-12 (A701-A712) primer tubes (orange caps, yellow solution) horizontally, aligned with columns 1 through 12.
- 5 Set out a new 96-well PCR plate (hereafter referred to as the **AMP** plate).
- 6 Add index primers to the AMP plate as follows:
 - a Add 9 μ l of the selected index primers [A–H (A501–A508)] to the appropriate well in a column of the **AMP** plate.
 - b Add 9 μ l of the selected index primers [1–12 (A701–A712)] to the appropriate row of the **AMP** plate.
Tips must be changed after each row to avoid index cross-contamination.

**NOTE**

Index primers are single use only and should be discarded immediately after use.

- 7 Prepare the PCR Master Mix/PCR Polymerase PCR working solution as follows:
 - a For 12 samples/24 libraries, add 21 μ l of PCR Polymerase to 1050 μ l of PCR Master Mix.
 - b Invert the prepared PCR working solution 20 times to mix.

The PCR working solution is stable at room temperature for 10 minutes and can also be stored on ice if needed during preparation.

Procedure

- 1 Remove the **FPU** from the incubator, and then remove the aluminum foil seal.
- 2 Cover the filter plate with the lid and centrifuge at $2400 \times g$ at room temperature for 2 minutes.
- 3 Add 25 μ l of 0.05 N NaOH to each sample well on the filter plate. Pipette the NaOH up and down 5–6 times to make sure the NaOH contacts the filter.

- 4 Cover and incubate the filter plate at room temperature for 5 minutes.
 - 5 While the filter plate is incubating, transfer 22 µl of the PCR working solution to each well of the AMP plate containing index primers.
 - 6 Transfer samples eluted from the filter to the AMP plate as follows:
 - a Pipette the samples in the first column of the filter plate up and down 5–6 times.
-  **NOTE**
Slightly tilt the **FPU** plate to ensure complete aspiration and to avoid air bubbles.
- b Transfer 20 µl from the filter plate to the corresponding column of the **AMP** plate.
 - c Gently pipette up and down 5–6 times to thoroughly combine the DNA with the PCR working solution.
 - d Transfer the remaining columns from the filter plate to the AMP plate in a similar manner. *Tips must be changed after each column to avoid index and sample cross-contamination.*
- 7 Seal the **AMP** plate and secure with a rubber roller.
 - 8 Centrifuge at 1000 × g at room temperature for 1 minute.
 - 9 Transfer the AMP plate to the post-amplification area.
 - 10 Perform PCR using the following program on a thermal cycler:
 - ▶ 95°C for 3 minutes
 - ▶ 32 cycles of:
 - ▶ 95°C for 30 seconds
 - ▶ 66°C for 30 seconds
 - ▶ 72°C for 60 seconds
 - ▶ 72°C for 5 minutes
 - ▶ Hold at 10°C

SAFE STOPPING POINT

If not proceeding immediately to PCR Clean-Up, the **AMP** plate can remain on the thermal cycler overnight, or can be stored at 2°C to 8°C up to 48 hours.

Verify Library Preparation

- 1 Combine 5 µl of amplified product with 15 µl of RNase/DNase-free water and DNA loading dye if necessary.
- 2 Run on a 2-4% TBE agarose gel with 100 bp ladder to confirm the presence and brightness of the 250-280 bp library product.
 - ▶ Samples that show amplification in both Oligo Pool A and B are considered valid and can be processed through the remainder of the workflow.
 - ▶ Samples showing little to no amplification in one Oligo Pool or both Oligo Pools are considered invalid and should not be processed through the remainder of the workflow. If an invalid gel result is observed, library preparation for that sample or samples will need to be repeated.
 - ▶ The Negative Control sample should show an Invalid gel result. If the Negative Control sample shows amplification in both Oligo Pool A and B, this indicates potential processing errors, i.e. contamination, pipetting errors, or sample mix-ups. The operator should evaluate the root cause and repeat library preparation as appropriate. If contamination is suspected, library preparation for all samples should be repeated. If pipetting errors or mix-ups are suspected, repeat library preparation for the affected samples.

NOTE

For more information and examples of valid and invalid results, see the *Praxis Extended RAS Panel Reference Guide* (document # 1000000037191).

PCR Clean-Up

Preparation

- 1 Bring the PCR Clean-Up Beads to room temperature.

- 2 Prepare fresh 80% ethanol from absolute ethanol. For 12 samples/24 libraries, prepare 80% ethanol using 12 mL of absolute ethanol and 3 mL of DNase/RNase-free water. Mix thoroughly.

Procedure

- 1 Centrifuge the **AMP** plate at $1000 \times g$ at room temperature for 1 minute.
- 2 Set out a new MIDI plate (hereafter referred to as the **CLP** plate).
- 3 Invert PCR Clean-Up Beads 10 times. Vortex vigorously and then invert 10 more times. Visually inspect solution to ensure that beads are resuspended.
- 4 Add 55 μ l of PCR Clean-Up Beads to each well of the **CLP** plate.
- 5 Transfer the entire PCR product from the **AMP** plate to the **CLP** plate.
- 6 Seal the **CLP** plate and shake on a microplate shaker at 1800 rpm for 2 minutes.
- 7 Incubate at room temperature without shaking for 10 minutes.
- 8 Place the plate on a magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
- 9 With the **CLP** plate on the magnetic stand, carefully remove and discard the supernatant.
- 10 With the **CLP** plate on the magnetic stand, wash the beads as follows:
 - a Add 200 μ l of freshly prepared 80% ethanol to each sample well.
 - b Incubate the plate on the magnetic stand for a minimum of 30 seconds or until the supernatant is clear.
 - c Carefully remove and discard the supernatant.
- 11 Repeat the wash as described in the previous step.
- 12 Use a P20 multi-channel pipette set to 20 μ l to remove excess ethanol. Do not disturb or touch the beads.
- 13 Remove the **CLP** plate from the magnetic stand and air-dry the beads for 5 minutes.
- 14 Add 25 μ l of Elution Buffer to each sample.
- 15 Seal the **CLP** plate and shake on a microplate shaker at 1800 rpm for 5 minutes. After shaking, verify if samples were resuspended. If not, repeat this step.
- 16 Place the **CLP** plate on the magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
- 17 Set out a new MIDI plate (hereafter referred to as the **LNP** plate).
- 18 Transfer 20 μ l of the supernatant from the **CLP** plate to the **LNP** plate.
- 19 [Optional] Transfer the remaining 5 μ l of supernatant from the **CLP** plate to a new plate and label the plate with a run name and date. Store this plate at -25°C to -15°C until completion of the sequencing run and data analysis. The cleaned up PCR products can be used for troubleshooting efforts in the event of sample failures.
- 20 Seal the **LNP** plate with an adhesive plate seal, and then centrifuge at $1000 \times g$ at room temperature for 1 minute to ensure all the supernatant is at the bottom of the well.

SAFE STOPPING POINT

After PCR Clean-Up, the plate is stable for up to 3 hours at 2°C to 8°C .

Library Normalization

Preparation

- 1 Prepare fresh 0.1 N NaOH by adding 30 μ l of 10 N NaOH to 2970 μ l of RNase/DNase-free water.
- 2 Bring Library Normalization Diluent, Library Beads, and Library Normalization Wash to room temperature.
- 3 Vortex Library Normalization Diluent vigorously and ensure that all precipitates have dissolved.
- 4 Vortex Library Beads vigorously for 1 minute with intermittent inversion until the beads are resuspended and no pellet is found at the bottom of the tube when the tube is inverted.

Procedure

- 1 Mix Library Normalization Diluent and Library Beads in a fresh 15 ml conical tube as follows:
 - a For 12 samples/24 libraries, add 1575 μ l of Library Normalization Diluent to a 15 ml conical tube.
 - b To homogenize beads, use a P1000 pipette set to 1000 μ l to resuspend Library Beads thoroughly by pipetting up and down 10 times.

**NOTE**

It is critical to completely resuspend the library bead pellet at the bottom of the tube. The use of a P1000 ensures that the beads are homogeneously resuspended and that there is no bead mass at the bottom of the tube. This is essential for achieving consistent cluster density on the flow cell.

- c For 12 samples/24 libraries, pipette 300 μ l of Library Beads to the tube containing Library Normalization Diluent.
 - d Mix by inverting the tube 15–20 times.
- 2 Add 45 μ l of the combined Library Normalization Diluent/Library Beads working solution to each well of the LNP plate containing libraries.
 - 3 Seal the LNP plate and shake on a microplate shaker at 1,800 rpm for 30 minutes.
 - 4 Place the LNP plate on a magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
 - 5 With the LNP plate on the magnetic stand, carefully remove and discard the supernatant.
 - 6 Remove the LNP plate from the magnetic stand and wash the beads with Library Normalization Wash as follows:
 - a Add 45 μ l of Library Normalization Wash to each sample well.
 - b Seal the LNP plate and shake on a microplate shaker at 1,800 rpm for 5 minutes.
 - c Place the LNP plate on the magnetic stand for a minimum of 2 minutes or until the supernatant is clear.
 - d Carefully remove and discard the supernatant.
 - 7 Repeat the Library Normalization Wash procedure as described in the previous step.
 - 8 Seal the LNP plate with an adhesive plate seal.
 - 9 Centrifuge the LNP plate at $1000 \times g$ at room temperature for 30 seconds to collect residual wash buffer.
 - 10 Place the LNP plate on the magnetic stand.
 - 11 Use a P20 multi-channel pipette set to 20 μ l to remove excess Library Normalization Wash.
 - 12 Remove the LNP plate from the magnetic stand and add 30 μ l of 0.1 N NaOH to each well.
 - 13 Seal the LNP plate and shake on a microplate shaker at 1,800 rpm for 5 minutes.
 - 14 During the 5-minute shake, set out a new 96-well PCR plate (hereafter referred to as the SGP plate.)
 - 15 Add 30 μ l Library Storage Buffer to each well to be used in the SGP plate.
 - 16 After the 5-minute shake, ensure that all samples in the LNP plate are resuspended. If the samples are not resuspended, gently pipette those samples up and down or lightly tap the plate on the bench to resuspend the beads, then shake for another 5 minutes.
 - 17 [Optional] Centrifuge the LNP plate at $1000 \times g$ at room temperature for 30 seconds to collect contents to the bottom of the wells.
 - 18 Place the LNP plate on the magnetic stand for a minimum of 2 minutes.
 - 19 Transfer the supernatant from the LNP plate to the SGP plate. Gently pipette up and down 5 times to mix.
 - 20 Seal the SGP plate and then centrifuge at $1,000 \times g$ at room temperature for 1 minute. Proceed immediately to Library Pooling.

Library Pooling

Prepare for Library Pooling

- 1 Set a heat block suitable for 1.5 ml centrifuge tubes to 96°C.
- 2 In an ice bucket, prepare an ice-water bath. Chill the Library Dilution Buffer in the ice-water bath.
- 3 Begin thawing the MiSeq™ Dx SBS cartridge v3.

Prepare the Reagent Cartridge

- 1 Thaw the MiSeq™ Dx SBS Cartridge v3 in a water bath containing enough room temperature deionized water to submerge the base of the reagent cartridge up to the water line printed on the reagent cartridge. Do not allow the water to exceed the maximum water line.
- 2 Allow the reagent cartridge to thaw in the room temperature water bath for approximately 1 hour or until thawed.

- 3 Remove the cartridge from the water bath and gently tap it on the bench to dislodge water from the base of the cartridge. Dry the base of the cartridge. Make sure that no water has splashed on the top of the reagent cartridge.

Denature and Dilute PhiX Internal Control

- 1 Prepare 0.1 N NaOH by combining the following volumes in a conical tube:

- ▶ DNase, RNase-free water (2475 µl)
- ▶ Stock 10 N NaOH (25 µl)

- 2 Invert the tube several times to mix.

CAUTION

Using freshly diluted NaOH is essential in order to denature PhiX for cluster generation on the MiSeqDx.

- 3 Combine the following volumes to make Diluted PhiX Internal Control:

- ▶ PhiX Internal Control (2 µl)
- ▶ 1X TE Buffer (8 µl)

- 4 Combine the following volumes to make Denatured Diluted PhiX Internal Control:

- ▶ Diluted PhiX Internal Control (10 µl)
- ▶ 0.1 N NaOH (10 µl)

- 5 Vortex briefly to mix the Denatured Diluted PhiX Internal Control.

- 6 Briefly centrifuge the Denatured Diluted PhiX Internal Control to collect contents.

- 7 Incubate for 5 minutes at room temperature to denature the PhiX Internal Control into single strands.

- 8 Combine the following volume of pre-chilled Library Dilution Buffer and Denatured Diluted PhiX Internal Control to make PhiX/Library Dilution Buffer solution.

- ▶ Denatured Diluted PhiX Internal Control (20 µl)
- ▶ Pre-chilled Library Dilution Buffer (980 µl)

Inspect the Reagent Cartridge

- 1 Invert the reagent cartridge ten times to mix the thawed reagents, and then inspect that all positions are thawed.

NOTE

It is critical that the reagents in the cartridge are thoroughly thawed and mixed to ensure proper sequencing.

- 2 Inspect reagents in positions 1, 2, and 4 to make sure that they are fully mixed and free of precipitates.

- 3 Gently tap the cartridge on the bench to reduce air bubbles in the reagents.

NOTE

The MiSeqDx sipper tubes go to the bottom of each reservoir to aspirate the reagents, so it is important that the reservoirs are free of air bubbles.

- 4 Place the reagent cartridge on ice or set aside at 2°C to 8°C (up to 6 hours) until ready to set up the run. For best results, proceed directly to loading the sample and setting up the run.

Prepare Samples for Sequencing

- 1 Bring Library Dilution Buffer to room temperature. Vortex Library Dilution Buffer and ensure that all precipitates have dissolved completely.

- 2 Centrifuge the SGP plate at 1000 × g at room temperature for 1 minute.

- 3 Set out a fresh microcentrifuge tube (hereafter referred to as the PAL tube.)

- 4 Determine the samples to be pooled for sequencing. A maximum of 12 samples/24 libraries can be pooled for each sequencing run.

- 5 Transfer 10 µl of each library to be sequenced from the SGP plate, column by column, to a PCR eight-tube strip. Seal SGP with an adhesive plate seal and set aside.

- 6 Combine and transfer the contents of the PCR eight-tube strip into the PAL tube. Mix the PAL tube thoroughly.

- 7 Set out 3–6 fresh microcentrifuge tubes (hereafter referred to as the DAL tube.)

- 8 Add 784 µl of Library Dilution Buffer to each DAL tube.

- 9 Add 8 µl of PhiX/Library Dilution Buffer solution to each **DAL** tube. Pipette up and down 3–5 times to rinse the tip and ensure complete transfer.
- 10 Transfer 8 µl of **PAL** to each **DAL** tube containing Library Dilution Buffer. Pipette up and down 3–5 times to rinse the tip and ensure complete transfer.
- 11 Mix the **DAL** tubes by vortexing the tube at top speed.

**NOTE**

The diluted library **DAL** can be stored at -25°C to -15°C for up to 30 days.

- 12 Centrifuge the one **DAL** tube at room temperature for 1 minute to collect contents.
- 13 Incubate the **DAL** tube on a heat block at 96°C for 2 minutes.

**WARNING**

Use caution when removing the tube. The tube and the heat block may be hot.

- 14 After the incubation, invert the **DAL** tube 1–2 times to mix, then immediately place in the ice-water bath.
- 15 Keep the **DAL** tube in the ice-water bath for 5 minutes.

**NOTE**

You must perform the heat denaturation step immediately before loading **DAL** into the MiSeq™ Dx SBS cartridge v3 ensure efficient template loading on the MiSeq™ Dx flow cell Replicate **DAL** tubes must be heat denatured immediately before MiSeqDx loading.

Library Sequencing

Procedure

For details on the steps outlined here, see the *MiSeqDx Instrument Reference Guide* (document # 1000000021961).

- 1 Use a separate, clean, and empty 1 ml pipette tip to pierce the foil seal over the reservoir on the MiSeqDx Cartridge labeled **Load Samples**.
- 2 Pipette 600 µl of the sample DAL libraries into the **Load Samples** reservoir. Take care to avoid touching the foil seal while dispensing the sample.
- 3 Check for air bubbles in the reservoir after loading sample. If air bubbles are present, gently tap the cartridge on the bench to release the bubbles.
- 4 Log in to the MiSeq Operating Software (MOS) using your Local Run Manager user name and password.

**NOTE**

If a previous user is logged in, select **Log Off** before logging in. Confirm that the proper username or login credentials are displayed.

- 5 From the Home screen, select **Sequence**.
- 6 Select a run from the list.
- 7 Confirm the run settings.
- 8 Clean the flow cell.
- 9 Load the flow cell.
- 10 Empty the waste bottle and load the MiSeq™ Dx SBS Solution (PR2) bottle.
- 11 Load the reagent cartridge.
- 12 The run automatically starts after passing pre-run checks.
- 13 Perform a post-run wash.

View Analysis Results

- 1 From the Local Run Manager dashboard, click the run name.
- 2 Click the Samples and Results tab to view the analysis report.
 - If analysis was repeated, expand the Select Analysis drop-down list and select the appropriate analysis.
- 3 [Optional] Click the **Copy to Clipboard**  icon to copy the Analysis folder path.

Interpretation of Results

- 1 The Illumina Praxis Extended RAS Panel is designed to detect the *KRAS* and *NRAS* mutations listed in Table 1.
- 2 Analysis results are summarized on the Samples and Results tab. The report is also available for export.

Praxis Extended RAS Panel Report

The Praxis Extended RAS Panel Report details if a run passed or failed, and whether each sample contains any of the 56 Extended RAS Panel mutations.

The Run State on the Samples and Results tab shows if the sequencing run passed quality metrics. If the Run State is Fail, all samples are listed as invalid.

Sample Results

Table 14 Sample Results Table

Column Heading	Description
Sample Name	The sample name provided when the run was created.
Sample	Sample validity results. Valid — Presence or absence of panel mutations is reported. Invalid — Presence or absence of panel mutations cannot be determined. See <i>Sample Invalid</i> on page 25.
Result	Sample results for panel mutation detection. See <i>Panel Mutation Detected</i> on page 25 and <i>Panel Mutation Not Detected</i> on page 25.
Gene	The gene, <i>KRAS</i> or <i>NRAS</i> , where a mutation is detected, if applicable.
Exon	The exon number in which the mutation is present, if applicable.
Amino Acid	Human Genome Variation Society (HGVS) amino acid change notation.
Nucleotide	HGVS nucleotide change notation.

Sample Invalid

The sample result, Invalid, indicates that the presence or absence of a Praxis Extended RAS Panel mutation cannot be determined. Repeat library preparation, sequencing, and analysis for Invalid samples. If a sample fails 2 times, additional troubleshooting can be necessary. If all samples in a run are Invalid because of a Run State of Fail, contact Illumina Technical Support.

Panel Mutation Detected

The sample result, Panel Mutation Detected, indicates that a Praxis Extended RAS Panel mutation is detected. The report lists all detected mutations.

Panel Mutation Not Detected

The sample result, Panel Mutation Not Detected, indicates that no Praxis Extended RAS Panel mutations are detected.

Quality Control Procedures

Good laboratory practices dictate that control material should be evaluated to detect differences in processing and technical procedures in the user's laboratory that may produce significant variability in results.

- 1 **Positive Controls (Panel Control)**—The Panel Control is required on every run. The Panel Control has 2 mutations, *KRAS* Exon 2 p.Gly12Asp c.35G>A and *NRAS* Exon 3 p.Gln61Lys c.181C>A. The Panel Control must generate the expected genotype to be valid. If the control is invalid, processing errors may have occurred. The software will fail the entire sequencing run and all samples will appear as Invalid.

- 2 **Negative (No Template/No DNA) Control**— The use of a negative (no template/no DNA) control (Blank) is required on every run. For a valid run, the Blank sample should appear as Invalid. If an unexpected outcome occurs in the Blank sample, the software will fail the entire sequencing run and all samples will appear as Invalid.
- 3 Before using this assay in the laboratory, the performance of the assay should be verified by testing a number of positive and negative samples with known performance characteristics.
- 4 All quality control requirements should be performed in accordance with local, state, and/or federal regulations or accreditation requirements.
- 5 If a sequencing run has failed, contact Illumina Technical Support to troubleshoot your run.

Performance Characteristics

Definitions Of Calculations Used In Performance Characteristics

Sample Qualification (Sample QC) Pass Rate: The Sample QC pass rate is defined as the percent of specimens that meet the sample quality control specification ($dCq \leq 5$).

Sample First Pass Rate: Sample first pass rate is defined as the percent of specimens that give a valid result the first time they are processed.

Agreement Calculations: The Praxis Extended RAS Panel is designed to specifically detect 56 mutations in the KRAS and NRAS genes. The assay simultaneously determines the presence or absence of each of the 56 mutations in every sequenced sample. The reference for agreement calculations is bi-directional Sanger sequencing or a control condition where indicated.

- 1 Overall percent agreement (OPA) is calculated as the number of mutant and non-mutant positions that agree with the reference divided by the total number possible mutant and non-mutant positions for all valid replicates. OPA is calculated on a per mutation basis; there are 56 mutations per sample. Overall percent agreement was calculated as follows:

$$(\# \text{ specimens} \times 56 \text{ mutant and non-mutant positions concordant with reference}) / (\# \text{ specimens} \times 56 \text{ mutant and non-mutant positions total}) = \text{OPA}$$
- 2 Positive Percent Agreement (PPA) is calculated as 100% times the number of true positives (mutant positions concordant with reference) divided by the sum of the number of true positives and the number of false negatives on a per mutation basis. Positive Percent Agreement was calculated as follows:

$$\# \text{ mutant positions} \times \# \text{ replicates} = \# \text{ potential mutant positions}$$

$$(\# \text{ true positives}) / (\# \text{ true positives} + \# \text{ false negatives}) = \text{PPA}$$
- 3 Negative Percent Agreement (NPA) is calculated as 100% times the number of true negatives (non-mutant positions concordant with reference) divided by the sum of the number of true negatives and the number of false positives. NPA is calculated on a per mutation basis; there are 56 mutations per sample. Negative Percent Agreement was calculated as follows:

$$(\# \text{ wild type specimens} \times \# \text{ replicates} \times 56 \text{ non-mutant positions}) + (\# \text{ mutant specimens} \times \# \text{ replicates} \times \# \text{ non-mutant positions}) = \# \text{ potential non-mutant positions}$$

$$(\# \text{ true negatives}) / (\# \text{ true negatives} + \# \text{ false positives}) = \text{NPA}$$
- 4 Percent Correct Call (PCC) is calculated as 100% times the number of correct calls divided by the number of calls attempted. This term is used for Reproducibility. The same calculation is described as percent agreement (positive or negative percent agreement) in other analytical studies.

Confidence Interval Calculations: One-sided 95% lower confidence bounds and 2-sided 95% confidence intervals were calculated using the Clopper-Pearson (exact binomial) method for analytical or clinical studies, or the Wilson Score method (where indicated) in clinical studies.

Analytical Accuracy – Comparison to an Orthogonal Method

Banked formalin-fixed, paraffin embedded (FFPE) colorectal cancer tissue samples collected during the Amgen panitumumab study (protocol no. 20050203 [T203 study]) were tested with the Praxis Extended RAS Panel. Samples were also evaluated using the reference method, bidirectional Sanger sequencing (Sanger). The agreements between results from Praxis Extended RAS Panel and Sanger were compared to evaluate the accuracy of the Praxis Extended RAS Panel.

Subject-level RAS status was positive if the subject had at least 1 mutation detected. Subject-level RAS status was negative if no mutation was detected. Positive Percent Agreement (PPA) is calculated as 100% times the number of true positive subjects divided by the sum of the number of true positive subjects and the number of false negative subjects. Negative Percent Agreement (NPA) is calculated as 100% times the number of true negative subjects divided by the sum of the number of true negative subjects and the number of false positive subjects. Two-sided 95% confidence intervals are calculated using the Clopper-Pearson (Exact binomial) method.

Of the 1183 subjects, the following were excluded from the accuracy analysis (combination of Sanger and Praxis Extended RAS Panel): 358 had no tissue available, 127 had insufficient tissue area, and 13 had insufficient neoplastic cells. Of the remaining 685 subjects, 214 were removed due to dCq > 5.0, 1 due to unestimable dCq, another 28 due to MiSeqDx QC failures, and lastly 1 subject was removed due to sample processing error/aborted run. The remaining 441 subjects had valid Sanger and Praxis Extended RAS Panel results, which were compared.

When evaluated at the subject level (Table 15), 230 of the 441 subjects were RAS positive by Sanger (mutation detected by Sanger); of these, 227 were RAS positive by the Praxis Extended RAS Panel. The remaining 211 of 441 subjects were RAS negative by Sanger (no mutation detected by Sanger); of these, 206 were RAS negative by the Praxis Extended RAS Panel. This resulted in a PPA of 98.7% and an NPA of 97.6%.

Table 15 Positive and Negative Percent Agreement of Subject-level RAS Results

Praxis Extended RAS Panel	Sanger		Total
	Positive	Negative	
Positive	227*	5	232
Negative	3**	206	209
Total	230	211	441

Performance Summary

Agreement Statistic	Point Estimate	Exact 95% CI
PPA	227/230 = 98.7%	[96.2%,99.7%]
NPA	206/211 = 97.6%	[94.6%,99.2%]

*There were 224 exact matches for within-subject, all mutation-level results; for 2 subjects, MiSeqDx detected the Sanger detected mutation and 1 additional mutation; for 1 subject, MiSeqDx and Sanger detected different mutations.

**1 subject had 2 mutations detected by Sanger; 2 subjects had 1 mutation detected by Sanger.

The data were also evaluated at the mutation level and designated true positive (TP), true negative (TN), false positive (FP), and false negative (FN). The PPA and NPA are presented for each mutation along with a 95% Clopper-Pearson (Exact binomial) confidence interval in Table 16 and Table 17.

Table 10 Positive and Negative Percent Agreement of Mutation-level Results in the KRAS Gene (dCq ≤ 5.0)

Location	Mutation	TP	FN	FP	TN	Total	PPA (95% CI)	NPA (95% CI)
KRAS Exon 2	c.34G>T	14	0	0	426	440	100.0 (76.8-100.0)	100.0 (99.1-100.0)
	c.34G>A	9	1	0	430	440	90.0 (55.5-99.7)	100.0 (99.1-100.0)
	c.34G>C	0	0	0	440	440	NE ¹	100.0 (99.2-100.0)
	c.34_35GG>TT	1	0	0	439	440	100.0 (2.5-100.0)	100.0 (99.2-100.0)
	c.34_35GG>AA	0	0	0	440	440	NE	100.0 (99.2-100.0)
	c.34_36GGT>TGG	0	0	0	440	440	NE	100.0 (99.2-100.0)
	c.35G>A	55	0	1	384	440	100.0 (93.5-100.0)	99.7 (98.6-100.0)
	c.35G>T	40	2	2	396	440	95.2 (83.8-99.4)	99.5 (98.2-99.9)
	c.35G>C	11	0	2	427	440	100.0 (71.5-100.0)	99.5 (98.3-99.9)
	c.37G>T	2	0	0	438	440	100.0 (15.8-100.0)	100.0 (99.2-100.0)
	c.37G>C	0	0	0	440	440	NE	100.0 (99.2-100.0)
	c.38G>A	41	0	0	399	440	100.0 (91.4-100.0)	100.0 (99.1-100.0)
	c.38_39GC>TT	0	0	0	440	440	NE	100.0 (99.2-100.0)
	c.38_39GC>AA	0	0	0	440	440	NE	100.0 (99.2-100.0)
	c.38_39GC>AT	0	0	0	440	440	NE	100.0 (99.2-100.0)
KRAS Exon 3	c.175G>A	1	0	0	440	441	100.0 (2.5-100.0)	100.0 (99.2-100.0)
	c.176C>G	1	0	0	440	441	100.0 (2.5-100.0)	100.0 (99.2-100.0)
	c.181C>A	0	1	0	440	441	0.0 (0.0-97.5)	100.0 (99.2-100.0)
	c.181C>G	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.182A>T	1	0	0	440	441	100.0 (2.5-100.0)	100.0 (99.2-100.0)
	c.182A>G	2	0	0	439	441	100.0 (15.8-100.0)	100.0 (99.2-100.0)
	c.183A>C	6	0	0	435	441	100.0 (54.1-100.0)	100.0 (99.2-100.0)
	c.183A>T	1	0	0	440	441	100.0 (2.5-100.0)	100.0 (99.2-100.0)
KRAS Exon 4	c.351A>C	2	1	0	437	440	66.7 (9.4-99.2)	100.0 (99.2-100.0)
	c.351A>T	2	0	0	438	440	100.0 (15.8-100.0)	100.0 (99.2-100.0)
	c.436G>A	11	0	0	429	440	100.0 (71.5-100.0)	100.0 (99.1-100.0)
	c.436G>C	1	0	0	439	440	100.0 (2.5-100.0)	100.0 (99.2-100.0)
	c.437C>T	3	0	0	437	440	100.0 (29.2-100.0)	100.0 (99.2-100.0)

* NE = not estimable

Table 17 Positive and Negative Percent Agreement of Mutation-level Results in the NRAS Gene (dCq ≤ 5.0)

Location	Mutation	TP	FN	FP	TN	Total	PPA (95% CI)	NPA (95% CI)
NRAS Exon 2	c.34G>T	1	0	1	439	441	100.0 (2.5-100.0)	99.8 (98.7-100.0)
	c.34G>A	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.34G>C	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.34_35GG>TT	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.34_35GG>AA	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.34_36GGT>TGG	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.35G>A	4	0	0	437	441	100.0 (39.8-100.0)	100.0 (99.2-100.0)
	c.35G>T	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.35G>C	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.37G>T	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.37G>C	2	0	0	439	441	100.0 (15.8-100.0)	100.0 (99.2-100.0)
	c.38G>A	0	0	1	440	441	NE	99.8 (98.7-100.0)
	c.38G>T	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.38_39GT>AA	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.38_39GT>AG	0	0	0	441	441	NE	100.0 (99.2-100.0)
NRAS Exon 3	c.175G>A	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.176C>G	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.181C>A	7	0	1	433	441	100.0 (59.0-100.0)	99.8 (98.7-100.0)
	c.181C>G	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.182A>T	2	0	0	439	441	100.0 (15.8-100.0)	100.0 (99.2-100.0)
	c.182A>G	5	0	0	436	441	100.0 (47.8-100.0)	100.0 (99.2-100.0)
	c.183A>C	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.183A>T	2	0	0	439	441	100.0 (15.8-100.0)	100.0 (99.2-100.0)
NRAS Exon 4	c.351G>C	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.351G>T	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.436G>A	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.436G>C	0	0	0	441	441	NE	100.0 (99.2-100.0)
	c.437C>T	0	0	0	441	441	NE	100.0 (99.2-100.0)

* NE = not estimable

Overall PPA and NPA calculated at the mutation level for each gene, and across both genes (total), are presented in Table 18. Two-sided 95% confidence intervals for PPA and NPA are calculated using the Clopper-Pearson (Exact binomial) method.


Table 18 Total Positive and Negative Percent Agreement of Mutation-Level Results

Location/ Mutation	TP	FN	FP	TN	Total	PPA (95% CI)	NPA (95% CI)
KRAS/All	204	5	5	12114	12328	97.6 (94.5-99.2)	99.96 (99.9-100.0)
NRAS/All	23	0	3	12322	12348	100 (85.2-100.0)	99.98 (99.9-100.0)
Total	227	5	8	24436	24676	97.8 (95.0-99.3)	99.97 (99.9-100.0)

Figure 7 shows a distribution of the Praxis Extended RAS Panel mutation frequencies. Using ≥ 50% tumor content in the FFPE specimen, the distribution of mutation frequencies observed in the '203 study for samples that tested positive is demonstrated in Figure 7. The majority of specimens showed mutation

frequencies in the 0.15 to 0.35 range. There were 14.9% of subjects with mutation frequency < 0.15, and 55.3% of subjects with mutation frequency ≥ 0.15 and ≤ 0.35 .

Praxis Extended RAS Panel Mutation Frequency Distribution for Positive Samples (n=235)

A cross-tabulation of Praxis Extended RAS Panel and Sanger mutation-level results are summarized in the confusion table of Sanger and Praxis Extended RAS Panel Results (see Table 19). There were 445 possible comparisons between mutations reported by Sanger and/or Praxis Extended RAS Panel for the 441 evaluable subjects in the RAS Evaluable analysis set. Four (4) subjects had two mutations reported by Praxis Extended RAS Panel, by Sanger, or both methods; therefore, the number of mutation-level results (n=445) is greater than the number of subjects (n=441).

Counts along the diagonal in the table represent results that agree at the mutation level (both methods report the same mutation or both methods report no mutation detected). Counts off the diagonal represent results that disagree. Both methods agreed on most results (434/445, 97.5%); 206 of the results that agreed were reported as no mutation detected (negative). In the negative column for Sanger, there were 5 instances (FP) of mutations reported by the Praxis Extended RAS Panel (eg, 1 sample Sanger reported as no mutation, while the Praxis Extended RAS Panel reported KRAS Ex2 c35C>A). In the negative row for Praxis Extended RAS Panel, there were 4 instances (FN) of mutations reported for Sanger (eg, 1 sample RAS Panel reported as no mutation, while Sanger reported KRAS Ex2 c34C>A).

Results corresponding to the 4 subjects who had 2 mutations detected by Sanger and/or Praxis Extended RAS Panel are noted by footnotes in the confusion table shown in Table 19.

Confusion Table of Sanger and Praxis Extended RAS Panel Results

		SANGER																					Negative	Total						
		KRAS Exon 2					KRAS Exon 3					KRAS Exon 4					NRAS Exon 2					NRAS Exon 3								
Missequenced RAS Panel		c34G>T	c34G>A	c34,35GG>TT	c35G>A	c35G>T	c35G>C	c37G>T	c38G>A	c178G>A	c178C>G	c181C>A	c182A>T	c182A>G	c183A>C	c183A>T	c351A>C	c351A>T	c436G>A	c436G>C	c437C>T	c34G>T	c35G>A	c37G>C	c38G>A	c181C>A	c182A>T	c182A>G	c183A>T	Negative
		14 ³	9	1	55	40	11	22	41	1	0	1	2	6	1	1	1	2	2	11	1	3	1	4	2	0	7	2	5	2
	c34G>T	14 ³																												14
	c34G>A		9																											9
	c34,35GG>TT			1																										1
	c35G>A				55																									55
	c35G>T					40																								40
	c35G>C						11																							11
	c37G>T							22																						22
	c38G>A								41																					41
	c178G>A									1																				1
	c178C>G										0																			0
	c181C>A											1																		1
	c182A>T												2																	2
	c182A>G													6																6
	c183A>C														1															1
	c183A>T															2														2
	c351A>C																1													1
	c351A>T																	2												2
	c436G>A																		11											11
	c436G>C																			1										1
	c437C>T																				3									3
	c34G>T																													1
	c35G>A																													4
	c37G>C																													2
	c38G>A																													1
	c181C>A																													8
	c182A>T																													1
	c182A>G																													5
	c183A>T																													2
	Negative																													206
	Total	15	10	1	55	42	11	22	42	1	0	1	2	6	1	1	3	2	11	1	3	7	5	2					211	445

Footnoted results A, B, C and D indicate the 4 subjects who contributed two results to the cell counts in the table

Praxis Extended RAS Panel

^ASubject where Sanger reported 2 mutations (*KRAS* Exon 2 c.34G>A and *KRAS* Exon 2 c.35G>T) and Praxis Extended RAS Panel reported panel mutation not detected.

^BSubject where Sanger reported 1 mutation (*KRAS* Exon 2 c.34G>T) and Praxis Extended RAS Panel reported the same mutation and an additional mutation (*KRAS* Exon 2 c.34G>T and *KRAS* Exon 2 c.35G>C).

^CSubject where Sanger reported 1 mutation (*KRAS* Exon 2 c.38G>A) and Praxis Extended RAS Panel reported the same mutation and an additional mutation (*KRAS* Exon 2 c.38G>A and *KRAS* Exon 2 c.35G>T).

^DSubject where Sanger reported two mutations (*KRAS* Exon 2 c.38G>A and *NRAS* Exon 2 c.35G>A) and Praxis Extended RAS Panel reported the same two mutations (*KRAS* Exon 2 c.38G>A and *NRAS* Exon 2 c.35G>A).

Analytical Sensitivity (DNA Input Range)

The Praxis Extended RAS Panel workflow involves sample qualification, library preparation and sequencing. Genomic DNA specimens are qualified for use in the assay by measuring amplifiability of the specimen relative to a control DNA template through a real-time quantitative PCR (qPCR) reaction. Measuring amplifiability accounts for both the quantity and quality of the DNA specimen for qualification. A quantification threshold value (Cq) is generated using the Quality Control Template and each genomic DNA specimen. The metric to assess the amplifiability is the change in quantification threshold (Cq) between the specimen and the control template DNA (Delta Cq; dCq). The dCq value is used to determine whether a specimen qualifies for the Praxis Extended RAS Panel assay and is the measure of DNA quality and quantity (DNA Input) for this assay. The acceptable measure for the assay is a Delta Cq ≤ 5 .

DNA input was evaluated over the range of delta Cq values from -0.5 to 5 in order to evaluate assay performance. A Delta Cq level of -0.5 to 5 corresponds to approximately 1200 ng to 25 ng, respectively, of intact DNA used in the assay.

A panel of 9 CRC FFPE specimens each containing 1 unique Extended RAS panel mutation, and 1 CRC FFPE specimen wild type for Extended RAS panel mutations was evaluated. The mutations represented were in *KRAS* exons 2, 3, and 4 and *NRAS* exons 2 and 3. The mutation frequency for this panel ranged from approximately 10% to 20%.

The specimen panel was tested in duplicate library preparations at each DNA input level for a total of 20 observations at each DNA input level. Specimens were sequenced at the maximum multiplexing level of 10 samples plus controls per run. Two (2) DNA input levels were tested: 1) dCq -0.5 to 0.5 (nominal high end) 2) dCq 4 to 5 (nominal low end).

All specimens showed a valid result after library preparation and sequencing with a sample first pass rate of 100% (20/20) for each DNA input level. OPA, PPA, and NPA are calculated on a per mutation basis; there are 56 mutations per sample. The overall percent agreement relative to Sanger sequencing was 100% (1120/1120) for each DNA input level. The positive percent agreement was 100% (18/18) and negative percent agreement was 100% (1102/1102) for both DNA input levels relative to Sanger sequencing. Therefore, this study demonstrates the performance of the Praxis Extended RAS Panel assay for the delta Cq (DNA input) range of -0.5 to 5.

To demonstrate the performance with challenging samples that harbor mutations near the limit of detection (LoD) of the mutant allele frequency, an additional study was conducted. Since DNA is qualified based on a dCq cutoff of 5.0, a nominal DNA input concentration range (12.5ng to 50ng) corresponding to approximate values of dCq 6-4 was selected for evaluation. A panel of 5 blended CRC FFPE specimen mixes, each containing 2 unique representative Extended RAS panel mutations, one at a low (approximately 8%) and one at a high (approximately 14%) mutation frequency level, was used. Two (2) of the five (5) sample mixes were spiked with a RAS panel multinucleotide variant (MNV). Two DNA input levels were evaluated for each FFPE blend: a dCq level of 4-5 corresponds to approximately 50-25 ng, respectively, and a dCq level of 5-6 corresponds to approximately 25-12.5 ng, respectively, of intact DNA. Twelve (12) observations of each of the 5 specimen mixes were collected over 3 non-consecutive days per DNA Input/Delta Cq and mutation frequency level. The total number of observations for the study across all specimen mixes was 60 observations; there were 120 observations across all mutation frequency levels per DNA input level (2 observations $\times$ 2 instruments [sequencing runs] $\times$ 3 non-consecutive days $\times$ 5 mixes). The sample first pass rate was 100% (60/60) for dCq 4 to 5 and 91.7% (55/60) for dCq 5 to 6. For dCq 4 to 5, 100% (120/120) of observations showed the expected Extended RAS Panel mutation result. For dCq 5 to 6, 95% (114/120) of observations showed the expected Extended RAS Panel mutation result. Mean mutation frequencies were similar between Delta Cq levels, whereas specimens with Delta Cq 5 to 6 had more mutation frequency variability in general than the same specimens with Delta Cq 4 to 5. The study confirmed that the Praxis Extended RAS Panel is able to detect mutants with accuracy at DNA input corresponding to dCq ≤ 5.0 and approximately 25 ng of intact DNA for mutations near the LoD.

Analytical Sensitivity (Limit of Blank [LoB] and Limit of Detection [LoD])

LoD is defined as the minor allele frequency at which 95% of samples are reliably detected. The Limit of Blank (LoB) and the Limit of Detection (LoD) study for the Praxis Extended RAS Panel was conducted using 3 blended CRC FFPE specimen mixes, each containing 2 representative Extended RAS Panel mutations at 5 mutation frequency levels (approx. 10%, 5%, 2.5%, 1.25% and 0%). These mixes were run in duplicate with 2

different reagent lots by 5 operators on 2 nonconsecutive days for a total of 20 observations per specimen mix per lot. The total number of observations for the study across all levels, specimen mixes, and reagent lots is 600 observations (20 observations × 5 levels × 3 specimen mixes × 2 lots = 600 observations). All runs included 10 samples plus controls on a single flow cell to test the LoD at the assay sample throughput limits. The sample mixes were prepared at the low level of DNA input (dCq=4–5). The range of results for 95% detection of each mutation tested is shown in Table 20. Using Probit analysis, the LoD for the 6 representative mutations ranged between 4.2%–5.3% (Table 20). Overall, the data supports the claim of 5%.

Table 20 Limit of Detection

Gene	Exon	Mutation	LoD
KRAS	2	c.35 G>A	5.0%
NRAS	3	c.182 A>G	4.7%
KRAS	2	c.38 G>A	4.8%
KRAS	4	c.436 G>A	5.3%
KRAS	3	c.183 A>C	4.2%
NRAS	2	c.38_39 GT>AA	4.5%

The limit of blank (LoB) for each of the 56 mutations claimed by this assay was calculated using specimens wild type for RAS panel mutations. Three wild type (WT) specimens were used in this study. A total of 60 observations were generated per study mutation per lot (20 observations × 3 WT samples). Two lots were used in the study. No false positive calls were reported. There was one observation with mutation frequency above the threshold, but the mutation was only detected in one oligo pool and not confirmed in the other, and was therefore not reported. The clinical threshold is set at 2.6%.

Precision

Lot-to-Lot Precision

A Lot-to-Lot Precision study was conducted to evaluate the performance of the Praxis Extended RAS Panel across manufactured reagent kit lots. Three (3) reagent kit lots were tested using a panel of 5 blended CRC FFPE specimens. Each CRC FFPE specimen contained 2 unique representative Extended RAS panel mutations: one at a low (approximately 8%) mutation frequency level and one at a high (approximately 14%) mutation frequency level. Twelve (12) observations of each of the 5 specimen mixes were collected over 3 nonconsecutive days with 3 reagent kit lots. The total number of observations for the study across all reagent lots was 180 observations across all specimen mixes and 360 observations across all mutation frequency levels. Across all lots and days, 99.7% (359/360) of observations showed the expected Extended RAS Panel mutation result. One low frequency mutation was incorrectly called as wild type. A variance component analysis was performed for each of the mutations/mutation frequency levels to estimate the variability of the system. The total standard deviations ranged from 0.011 to 0.029. The reagent lot component of the total standard deviation ranged from 0 to 0.015.

Reproducibility

An external study was performed to assess the reproducibility of the Praxis Extended RAS Panel across 3 external testing sites (2 operators per site), 1 reagent lot, and 3 non-consecutive testing days. Testing was conducted with 6 well-characterized sample panels of genomic DNA samples from FFPE clinical specimens or cell lines. Each panel consisted of 10 members, for a total of 60 members across panels.

The 60 panel members consisted of duplicates of 4 unique wild type (for RAS panel mutations) specimens, 12 unique mutant specimens (with a single mutation) prepared at both high and low mutation frequency levels, and 2 unique mutant specimens (with a single mutation) prepared at a low mutation frequency level only. Each unique specimen/mutation frequency level sample (tested in duplicate in each run) had 36 possible results (2 replicates × 2 operators × 3 days × 3 sites) if all results were valid. The Delta Cq values for each of the reproducibility sample panel members ranged from 3.64 to 4.08.

Percent Correct Call (PCC) across mutations and wild types was evaluated by comparing the Praxis Extended RAS Panel result to the expected mutation result (expected mutation detected or not detected) in each sample. PCC is calculated as 100% times the number of correct calls divided by the number of calls attempted. Two-sided 95% confidence intervals are calculated using the Wilson score method.

Qualitative Reproducibility

Combining sites, sample pass rates were $\geq 94.7\%$ for the first run of the sample or in samples tested in runs that were valid on first pass. The mutation-level percent correct call (PCC) across mutations from all mutant samples was 99.6% (905/909) excluding failed sample data. The number of attempted calls across all 56 panel mutations (regardless of whether a detected mutation was expected or not) for all valid samples was 58,856 (56×1051). Out of these 58,856 mutation level observations, there were only 6 incidents where the observed and expected results were discordant; the mutation-level PCC across mutations and wild types from all samples combined was 99.99% (58,850/58,856).

Table 21 shows PCCs for each of the 26 unique mutant specimen/mutation frequency level samples (combining replicates across 2 panel members) and corresponding 95% 2-sided score CIs.

Table 21 Mutation-Level Percent Correct Call for Individual Unique Mutant Specimens/Mutation Frequency Level Samples (Evaluable Analysis Set)

Mutation Level	Exon	Amino Acid Change	Base Change	N	Percent Correct Call (%)	95% CI ¹
Low	KRAS exon 2	p.Gly12Cys	c.34G>T	36	36/36(100%)	(90.4, 100)
		p.Gly12Arg	c.34G>C	36	36/36(100%)	(90.4, 100)
		p.Gly12Ser	c.34G>A	35	35/35(100%)	(90.1, 100)
		p.Gly12Ala	c.35G>C	35	35/35(100%)	(90.1, 100)
		p.Gly12Asp	c.35G>A	36	36/36(100%)	(90.4, 100)
		p.Gly12Val	c.35G>T	36	36/36(100%)	(90.4, 100)
		p.Gly13Asp	c.38G>A	35	35/35(100%)	(90.1, 100)
	KRAS exon 2(all samples)			249	249/249(100%)	(98.5, 100)
	KRAS exon 3	p.Gln61His	c.183A>C	33	33/33(100%)	(89.6, 100)
	KRAS exon 4	p.Lys117Asn	c.351A>T	35	33/35(94.3%)	(81.4, 98.4)
		p.Ala146Thr	c.436G>A	36	36/36(100%)	(90.4, 100)
	KRAS exon 4(all samples)			71	69/71(97.2%)	(90.3, 99.2)
	KRAS(all exons)			353	351/353(99.4%)	(98, 99.8)
	NRAS exon 2	p.Gly12Asp	c.35G>A	31	31/31(100%)	(89, 100)
	NRAS exon 3	p.Gln61Arg	c.182A>G	33	33/33(100%)	(89.6, 100)
	NRAS exon 4	p.Lys117Asn	c.351G>T	35	35/35(100%)	(90.1, 100)
		p.Ala146Thr	c.436G>A	35	34/35(97.1%)	(85.5, 99.5)
	NRAS exon 4(all samples)			70	69/70(98.6%)	(92.3, 99.7)
	NRAS(all exons)			134	133/134(99.3%)	(95.9, 99.9)
	RAS(all exons)			487	484/487(99.4%)	(98.2, 99.8)
High	KRAS exon 2	p.Gly12Arg	c.34G>C	33	33/33(100%)	(89.6, 100)

Mutation Level	Exon	Amino Acid Change	Base Change	N	Percent Correct Call (%)	95% CI ¹
		p.Gly12Ala	c.35G>C	36	36/36(100%)	(90.4, 100)
		p.Gly12Asp	c.35G>A	36	36/36(100%)	(90.4, 100)
		p.Gly12Val	c.35G>T	35	35/35(100%)	(90.1, 100)
		p.Gly13Asp	c.38G>A	36	36/36(100%)	(90.4, 100)
	KRAS exon 2(all samples)			176	176/176(100%)	(97.9, 100)
	KRAS exon 3	p.Gln61His	c.183A>C	35	35/35(100%)	(90.1, 100)
	KRAS exon 4	p.Lys117Asn	c.351A>T	36	36/36(100%)	(90.4, 100)
		p.Ala146Thr	c.436G>A	35	35/35(100%)	(90.1, 100)
	KRAS exon 4(all samples)			71	71/71(100%)	(94.9, 100)
	KRAS(all exons)			282	282/282(100%)	(98.7, 100)
	NRAS exon 2	p.Gly12Asp	c.35G>A	34	34/34(100%)	(89.8, 100)
	NRAS exon 3	p.Gln61Arg	c.182A>G	34	34/34(100%)	(89.8, 100)
	NRAS exon 4	p.Lys117Asn	c.351G>T	36	35/36(97.2%)	(85.8, 99.5)
		p.Ala146Thr	c.436G>A	36	36/36(100%)	(90.4, 100)
	NRAS exon 4(all samples)			72	71/72(98.6%)	(92.5, 99.8)
	NRAS(all exons)			140	139/140(99.3%)	(96.1, 99.9)
	RAS(all exons)			422	421/422(99.8%)	(98.7, 100)
Overall	All panel mutations			909	905/909(99.6%)	(98.9, 99.8)

^[1] 95% two-sided confidence interval calculated using the Wilson score method.

Table 20 Sample-Level Percent Correct Call Across Wild Type Samples (Evaluable Analysis Set)

Sample	N	Percent Correct Call (%)	95% CI ¹
Wild Type 1	35	35/35(100%)	(90.1, 100)
Wild Type 2	36	36/36(100%)	(90.4, 100)
Wild Type 3	36	36/36(100%)	(90.4, 100)
Wild Type 4	35	35/35(100%)	(90.1, 100)
Overall	142	142/142(100%)	(97.4, 100)

^[1] 95% two-sided confidence interval calculated using the Wilson score method.

Table 23 shows the variance component analysis of mutation frequency results for each source of variation and total variation for the Evaluable analysis set (for mutation frequency). Results are shown for each unique mutant specimen/mutation frequency level sample, combining sites.

Table 23: Variance Component Estimates for Individual Mutant Specimen and Mutation Level (Evaluable Analysis Set for Mutation Frequency)

Mutation Level	Exon	Amino Acid Change	Base Change	N	Mean Mutation Frequency	Site SD (%CV)	Operator (Site) SD (%CV)	Day(Site, Operator) SD (%CV)	Replicate SD (%CV)	Total SD (%CV)
Low	KRAS exon 2	p.Gly12Cys	c.34G>T	36	0.080	0.003 (4.0)	0.003 (4.0)	0.000 (0.0)	0.009 (10.8)	0.010 (12.2)
		p.Gly12Arg	c.34G>C	36	0.078	0.020 (25.2)	0.004 (5.3)	0.005 (6.6)	0.008 (9.8)	0.022 (28.3)
		p.Gly12Ser	c.34G>A	35	0.080	0.005 (6.1)	0.000 (0.0)	0.000 (0.0)	0.011 (13.9)	0.012 (15.2)
		p.Gly12Ala	c.35G>C	35	0.072	0.000 (0.0)	0.000 (0.0)	0.010 (13.8)	0.013 (17.5)	0.016 (22.2)
		p.Gly12Asp	c.35G>A	36	0.072	0.000 (0.0)	0.000 (0.0)	0.002 (3.4)	0.013 (18.7)	0.014 (19.0)
		p.Gly12Val	c.35G>T	36	0.083	0.002 (3.0)	0.000 (0.0)	0.000 (0.0)	0.009 (11.1)	0.010 (11.5)
		p.Gly13Asp	c.38G>A	35	0.092	0.000 (0.0)	0.000 (0.0)	0.015 (16.7)	0.017 (18.9)	0.023 (25.3)
	KRAS exon 3	p.Gln61His	c.183A>C	33	0.086	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.011 (12.9)	0.011 (12.9)
	KRAS exon 4	p.Lys117Asn	c.351A>T	33	0.090	0.000 (0.0)	0.000 (0.0)	0.003 (3.1)	0.010 (10.7)	0.010 (11.2)
		p.Ala146Thr	c.436G>A	36	0.099	0.005 (4.8)	0.000 (0.0)	0.000 (0.0)	0.009 (9.4)	0.010 (10.6)
	NRAS exon 2	p.Gly12Asp	c.35G>A	31	0.078	0.002 (2.9)	0.000 (0.0)	0.000 (0.0)	0.009 (12.1)	0.010 (12.4)
	NRAS exon 3	p.Gln61Arg	c.182A>G	33	0.097	0.000 (0.0)	0.000 (0.0)	0.002 (1.8)	0.009 (9.5)	0.009 (9.6)
	NRAS exon 4	p.Lys117Asn	c.351G>T	35	0.086	0.001 (1.7)	0.002 (2.4)	0.003 (3.9)	0.007 (8.4)	0.008 (9.7)
		p.Ala146Thr	c.436G>A	35	0.080	0.006 (6.9)	0.000 (0.0)	0.000 (0.0)	0.016 (20.5)	0.017 (21.7)
High	KRAS exon 2	p.Gly12Arg	c.34G>C	33	0.230	0.047 (20.4)	0.004 (1.7)	0.000 (0.0)	0.019 (8.1)	0.051 (22.1)

Mutation Level	Exon	Amino Acid Change	Base Change	N	Mean Mutation Frequency	Site SD (%CV)	Operator (Site) SD (%CV)	Day(Site, Operator) SD (%CV)	Replicate SD (%CV)	Total SD (%CV)
High	KRAS exon 2	p.Gly12Ala	c.35G>C	36	0.281	0.009 (3.3)	0.000 (0.0)	0.008 (2.9)	0.013 (4.7)	0.018 (6.5)
		p.Gly12Asp	c.35G>A	36	0.259	0.000 (0.0)	0.000 (0.0)	0.005 (2.0)	0.015 (5.8)	0.016 (6.2)
		p.Gly12Val	c.35G>T	35	0.206	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.020 (9.9)	0.020 (9.9)
		p.Gly13Asp	c.38G>A	36	0.263	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.024 (9.3)	0.024 (9.3)
	KRAS exon 3	p.Gln61His	c.183A>C	35	0.259	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.013 (5.0)	0.013 (5.0)
	KRAS exon 4	p.Lys117Asn	c.351A>T	36	0.267	0.008 (3.0)	0.007 (2.6)	0.003 (1.2)	0.014 (5.3)	0.018 (6.7)
		p.Ala146Thr	c.436G>A	35	0.277	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.015 (5.6)	0.015 (5.6)
	NRAS exon 2	p.Gly12Asp	c.35G>A	34	0.274	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.013 (4.8)	0.013 (4.8)
	NRAS exon 3	p.Gln61Arg	c.182A>G	34	0.281	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.012 (4.3)	0.012 (4.3)
	NRAS exon 4	p.Lys117Asn	c.351G>T	36	0.258	0.010 (3.9)	0.010 (3.7)	0.012 (4.5)	0.045 (17.3)	0.048 (18.7)
		p.Ala146Thr	c.436G>A	36	0.258	0.000 (0.0)	0.000 (0.0)	0.000 (0.0)	0.014 (5.6)	0.014 (5.6)

Analytical Specificity (Interfering Substances)

An Interfering Substances study was conducted to evaluate the impact of potentially interfering substances on the performance of the Praxis Extended RAS Panel Assay. Assay performance was evaluated in the presence of exogenous substances, (paraffin wax, xylene, ethanol and Proteinase K, extraction solutions), as well as endogenous substances (necrotic tissue and hemoglobin).

Exogenous Substances

The exogenous substances tested are extraction solutions commonly used during the DNA extraction process and are listed with tested quantities in Table 24. Fifteen (15) colorectal FFPE specimens were tested per interfering substance and compared with untreated controls. The specimens represented wild type specimens containing no Extended RAS panel mutations (5/15 specimens), as well as specimens containing prevalent Extended RAS panel mutations (10/15 specimens). Specimens were sequenced at the maximum multiplexing level of 10 samples plus controls per run.

Substances Tested

Interfering Substance	Actual Amount [ul/25uL eluate]
Deparaffinization Solution	1.69×10^{-04}
Paraffin Wax (in Xylene)	2.50×10^{-05}
Xylene	2.50×10^{-05}
Ethanol	1.69×10^{-04}
Proteinase K	3.30×10^{-06}
Wash Solution, N	6.25×10^{-01}
1X Wash Solution, P	6.25×10^{-01}
Wash Buffer AW1	6.25×10^{-02}
Wash Buffer AW2	6.25×10^{-01}

For all exogenous substances tested, all 15 specimens passed the sample qualification requirement (15/15, 100% Sample QC Pass Rate) and showed a valid result after library preparation and sequencing (15/15, 100% Sample First Pass Rate).

OPA, PPA, and NPA are calculated on a per mutation basis; there are 56 mutations per sample. All 15 specimens for all 9 exogenous substances showed agreement with the untreated control condition at all mutant and non-mutant positions. No false positives and no false negatives were observed. None of the potentially interfering substances evaluated at the maximum concentrations expected to be encountered in the process of genomic DNA (gDNA) extraction from FFPE tissue impacts performance of the Praxis Extended RAS Panel Assay.

Endogenous substances (Hemoglobin)

Fifteen (15) CRC FFPE specimens were each tested in the presence or absence of 2 mg/mL hemoglobin, a CLSI "high" amount of hemoglobin. The specimens represented wild type specimens containing no Extended RAS panel mutations (5/15 specimens), as well as specimens containing prevalent Extended RAS Panel mutations (10/15 specimens). Specimens were sequenced at the maximum multiplexing level of 10 samples plus controls per run. All 15 specimens passed the sample qualification requirement (15/15, 100% Sample QC Pass Rate) and showed a valid result after library preparation and sequencing (15/15, 100% Sample First Pass Rate). All 15 specimens showed agreement at all mutant and non-mutant positions with the untreated control condition. No false positives or false negatives were observed. The concentration of hemoglobin tested does not affect performance of the Praxis Extended RAS Panel.

Necrotic tissue

Fifteen (15) CRC FFPE specimens composed of wild type specimens containing no Extended RAS panel mutations (10/15 specimens), as well as specimens containing prevalent Extended RAS Panel mutations (5/15 specimens) and between 10–80% necrotic tissue as determined by pathology review, were utilized for evaluation of endogenous necrotic specimens. Specimens were sequenced at the maximum multiplexing level of 10 samples plus controls per

Praxis Extended RAS Panel

run. 14/15 specimens gave a valid result after library preparation and sequencing (93.3% Sample First Pass Rate). The overall percent agreement was 99.9% (783/784) relative to Sanger sequencing. PPA was 100% (4/4) and NPA was 99.9% (779/780). The 1 false positive detected was likely due to a mutation frequency below the limit of detection of Sanger sequencing. Overall, the Praxis Extended RAS Panel meets performance characteristics with tissue containing 10–80% necrosis.

Guard banding Studies

The Praxis Extended RAS Panel workflow guard banding studies evaluated selected workflow steps of Sample Qualification, Library Preparation, and Sequencing. Test samples included FFPE genomic DNA and cell line genomic DNA to increase mutation diversity. The Library preparation conditions were tested at the maximum multiplexing level of the assay at 10 specimens per run.

Sample Qualification (qPCR step)

Five genomic DNA samples were tested at various storage times of qPCR master mix and reaction mix. Since both the qPCR master mix and reaction mix are used right away, the testing conditions included holding the qPCR master mix at room temperature for 1 hour in ambient light, and the reaction mix for 10 minutes on ice. Average dCq differences of less than 0.25 were observed indicating stability of the qPCR master mix and reaction mixes.

Library Preparation

Hybridization and Extension-Ligation incubation times and temperatures were tested at approximately +/- 5% of the specified value. 100% of specimen tested showed a valid result with the expected mutations.

During PCR Amplification, sodium hydroxide is used to remove the extension-ligation product from the filter plate. The concentration and incubation times were tested at approximately +/- 5% of the specified value. 90% of specimens tested at low sodium hydroxide concentration and 100% of specimen tested at all other conditions showed a valid result with the expected mutations.

Storage times of PCR working solution and hold times of the PCR plate (post-amplification) were tested at approximately +/- 5% of the specified value. All specimens tested at all conditions showed a valid result with the expected mutations.

For the PCR clean-up and Library Normalization steps, the following conditions were tested at approximately +/- 5% of the specified value: a) incubation times, b) shaking times, and concentrations of c) ethanol (PCR clean-up), d) sodium hydroxide solution (Library Normalization). 90% of specimens tested for the longer PCR clean up elution incubation and 100% of specimens tested at all other conditions showed a valid result with the expected mutations.

For the Library Pooling step, times for heat denaturation and cooling were tested at approximately +/- 5% of the specified values. 90% of specimens tested at the longer cooling, and 100% of specimens tested at all other heat denaturation and cooling conditions showed a valid result with the expected mutations.

DAL (Diluted Amplicon Library) storage was tested up to 84 days (stored at -20°C). 100% of specimens tested showed a valid result with the expected mutations.

Sequencing

Robustness of deep sequencing is achieved with a minimum sequencing depth requirement for the samples to be called "valid". A total of 40 specimens plus controls were tested by 2 operators for meeting the minimum coverage of 1800X across both pools (900x per pool) per sample at the maximum multiplexing level of the assay at 10 specimens per run. All specimens met the minimum coverage requirement with an average coverage of approximately 31000x across both pools per sample.

Stability

Real-Time Stability

The objective of the real-time stability study was to establish the shelf life of the Praxis Extended RAS Panel kit by evaluating the product throughout the proposed shelf life when stored per label conditions. The study design was based on the testing of three Praxis Extended RAS Panel kit lots and utilized the classic stability study design described in CLSI EP25-A. The kits were stored in its final kitted configuration for the duration of the study, at storage conditions defined on the product label. Frozen kit components were stored at -15°C to -25°C. Refrigerated kit components were stored at 2°C to 8°C. Room-temperature components were stored at 15°C to 30°C. Each kit has the appropriate expiry date printed on the box level label.

Open-Tube Stability

The open-tube, in-use stability of the sample qualification and library preparation reagents was evaluated under standard use conditions. The reagents are expected to be used two (2) times which corresponds to opened tubes for the room temperature and 2°C to 8°C components and two (2) freeze/thaw cycles for the reagents stored between -15°C and -25°C. The frozen reagents were subjected to 5 and 6 freeze/thaw cycles between execution control testing and final testing. Room temperature and 2°C to 8°C components were tested under the standard use conditions per label requirements for two separate library preparations with opened reagents closed and returned to their indicated storage conditions between uses. The study showed that all quality control acceptance criteria were met for all kits tested, both in the execution tests and after the freeze/thaw cycles.

Summary of Analytical Studies

When evaluated at the specimen level at external sites using study specimens, and after exclusion of Praxis invalids, the Praxis Extended RAS Panel demonstrated a PPA of 98.7% (227/230) and an NPA of 97.6% (206/211), as compared to Sanger.

The Praxis Extended RAS Panel can detect 56 mutations () in the KRAS and NRAS genes. The performance of the Praxis Extended RAS Panel assay was demonstrated for the delta Cq (DNA input) range of -0.5 to 5. The LoD for each mutation tested ranged from 4.2% to 5.3%. Assay performance was not negatively impacted by interfering substances and necrotic tissue. No false positive results were detected due to well-to-well or run-to-run sample carryover. Precision and reproducibility of assay performance was demonstrated across multiple sources of variation with minimal contribution of reagent lot and other factors to overall variability.

References

- 1 Sorich MJ, Wiese MD, Rowland A, Kichenadasse G, McKinnon RA, Karapetis CS. Extended RAS mutations and anti-EGFR monoclonal antibody survival benefit in metastatic colorectal cancer: a meta-analysis of randomized, controlled trials. *Ann Oncol*. 2015 Jan;26(1):13-21.
- 2 Sepulveda AR, Hamilton SR, Allegra CJ, Grody W, Cushman-Vokoun AM, Funkhouser WK, et al. Molecular Biomarkers for the Evaluation of Colorectal Cancer: Guideline From the American Society for Clinical Pathology, College of American Pathologists, Association for Molecular Pathology, and American Society of Clinical Oncology. *J Mol Diagn* 2017;19:187-225.
- 3 Vaughn CP, Zobell SD, Furtado LV, Baker CL, Samowitz WS. Frequency of *KRAS*, *BRAF*, and *NRAS* mutations in colorectal cancer. *Genes Chromosomes Cancer*. 2011 May;50(5):307-12.
- 4 US NCCN Guidelines: Colon Cancer V.1.2017.

Patents and Trademarks

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2017 Illumina, Inc. All rights reserved.

Illumina, Genetic Energy, MiSeqDx, Powered by Illumina, Praxis, TruSeq, and the streaming bases design are registered or pending trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

This product may be used as a component of human molecular diagnostic assays where applicable country laws allow. This product by itself does not provide any diagnostic result. This product includes GoTaq® Hot Start Polymerase manufactured by Promega Corporation for distribution by Illumina, Inc. Licensed to Promega Corporation under U.S. Patent Nos. 5,338,671 and 5,587,287 and their corresponding foreign patents.

NOTICE TO PURCHASER: LIMITED LICENSE

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,994,056 and 6,171,785. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product solely in the field of quality assurance and quality control for products and processes in the In Vitro Human Diagnostics Field, including, without limitation, conformance with specifications, purity and batch-to-batch consistency and including reporting results of purchaser's activities for a fee or other commercial consideration. No right under any other patent claim is conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

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- 1 Soric MJ, Wiese MD, Rowland A, Kichenadasse G, McKinnon RA, Karapetis CS. Extended RAS mutations and anti-EGFR monoclonal antibody survival benefit in metastatic colorectal cancer: a meta-analysis of randomized, controlled trials. *Ann Oncol*. 2015 Jan; 26(1): 13-21.
- 2 Sepulveda AR, Hamilton SR, Allegra CJ, Grody W, Cushman-Vokoun AM, Funkhouser WK, et al. Molecular Biomarkers for the Evaluation of Colorectal Cancer: Guideline From the American Society for Clinical Pathology, College of American Pathologists, Association for Molecular Pathology, and American Society of Clinical Oncology. *J Mol Diagn* 2017; 19: 187-225.
- 3 Vaughn CP, Zobel SD, Furtado LV, Baker CL, Samowitz WS. Frequency of *KRAS*, *BRAF*, and *NRAS* mutations in colorectal cancer. *Genes Chromosomes Cancer*. 2011 May; 50(5): 307-12.
- 4 US NCCN Guidelines: Colon Cancer V.1.2017.

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
personally appeared Aurelia Brownridge,
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are
subscribed to the within instrument and acknowledged to me that he/she/they executed the same in
his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the
person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

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Local Run Manager Praxis Extended RAS Panel Analysis Module Workflow Guide

FOR IN VITRO DIAGNOSTIC USE

Overview	3
Set Parameters	3
View Analysis Results	4
Analysis Report	4
Analysis Output Files	5
Technical Assistance	6

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© 2017 Illumina, Inc. All rights reserved.

Illumina, Praxis, and the streaming bases design are registered or pending trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Overview

The Local Run Manager Praxis™ Extended RAS Panel analysis module is for use with Praxis Extended RAS Panel. The analysis module performs secondary analysis and report generation from sequencing runs that use the Praxis Extended RAS Panel. See the package insert *Praxis Extended RAS Panel* (document # 1000000037190).

About This Guide

This guide provides instructions for setting up run parameters for sequencing and analysis for the Praxis Extended RAS Panel analysis module. For information about the Local Run Manager dashboard and system settings, see the *Local Run Manager Software Reference Guide* (document # 1000000011880).

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Enter Samples Manually

- 1 Click the plate region A1 or E7 to enter sample information in the samples table.
Region A1 is equivalent to the upper left quadrant of wells in the Hyb plate layout, and includes samples for Oligo Pools A and B.
Region E7 is equivalent to the lower right quadrant of wells in the Hyb plate layout, and includes samples for Oligo Pools A and B.
- 2 In the Pool A table, enter a unique sample name in the Sample Name field for each well you use.
Use alphanumeric characters, dashes, or underscores.
The sample name automatically populates the corresponding well in the Pool B table.
Dual index adapters for each Index Read, specific to the well, display as you enter a sample name.
- 3 [Optional] Although the default layout is recommended, to change the position of control samples, use the Positive Control or Negative Control drop down menu.
- 4 Click the print  icon to display the plate layout.

- 5 Select **Print** to print the plate layout as a reference for preparing libraries.
- 6 Click **Close**.
- 7 Click **Save Run**.

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View Analysis Results

- 1 From the Local Run Manager dashboard, click the run name.
- 2 Click the **Samples and Results** tab to view the analysis report.
 - If analysis was repeated, expand the **Select Analysis** drop-down list and select the appropriate analysis.
- 3 [Optional] Click the **Copy to Clipboard**  icon to copy the Analysis folder path.

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The Run State on the **Samples and Results** tab shows if the sequencing run passed quality metrics. If the Run State is Fail, all samples are listed as invalid.

Sample Results

Table 1 Sample Results Table

Column Heading	Description
Sample Name	The sample name provided when the run was created.
Sample	Sample validity results. Valid – Presence or absence of panel mutations is reported. Invalid – Presence or absence of panel mutations cannot be determined. See <i>Sample Invalid</i> on page 5.
Result	Sample results for panel mutation detection. See <i>Panel Mutation Detected</i> on page 5 and <i>Panel Mutation Not Detected</i> on page 5.
Gene	The gene, <i>KRAS</i> or <i>NRAS</i> , where a mutation is detected, if applicable.
Exon	The exon number in which the mutation is present, if applicable.
Amino Acid	Human Genome Variation Society (HGVS) amino acid change notation.
Nucleotide	HGVS nucleotide change notation.

Sample Invalid

The sample result, **Invalid**, indicates that the presence or absence of a Praxis Extended RAS Panel mutation cannot be determined. Repeat library preparation, sequencing, and analysis for **Invalid** samples. If a sample fails 2 times, additional troubleshooting can be necessary. If all samples in a run are **Invalid** because of a Run State of Fail, contact Illumina Technical Support.

Panel Mutation Detected

The sample result, **Panel Mutation Detected**, indicates that a Praxis Extended RAS Panel mutation is detected. The report lists all detected mutations.

Panel Mutation Not Detected

The sample result, **Panel Mutation Not Detected**, indicates that no Praxis Extended RAS Panel mutations are detected.

Analysis Output

The following analysis output file is generated for the Praxis Extended RAS Panel analysis module. Analysis output files are located in the Alignment folder.

File Name	Description
ExtendedRasPanelReport.txt	The Praxis Extended RAS Panel report details if a run passed or failed, and whether each sample contains any of the 56 Extended RAS Panel mutations.

Alignment Folders

Each time that analysis is requested, the Local Run Manager creates an Alignment folder in the run folder named **Alignment_N**, where N is a sequential number.

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Reagent kit Praxis™ Extended RAS Panel for qualitative determination of the status of specific mutations genes of RAS (*KRAS* (exons 2, 3, and 4) and *NRAS* (exons 2, 3, and 4) by sequencing a new generation in a human genomic DNA sample isolated from colorectal cancer tissue samples fixed in formalin and enclosed in paraffin blocks

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

State of California
County of San Diego

On June 16, 2020 before me, Mary Agostino, Notary Public
(insert name and title of the officer)

personally appeared Angelina Brownridge
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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Компонент/ Номер по каталогу (PN)	Фотографии	Количество, объем/заполнени е, размеры: высота (H) длина (L) диаметр (D)	Активные ингредиенты, pH и физические и химические свойства
Инструкция по применению PN 15054121	 Praxis® Extended RAS Panel Praxis® Extended RAS Panel is a quantitative, colorimetric test using targeted high throughput parallel sequencing to detect the presence of specific mutations in RAS genes (KRAS, NRAS, HRAS, and BRAF). The test is performed on a 96-well plate and is suitable for use in a laboratory setting. The Praxis® Extended RAS Panel is a quantitative, colorimetric test using targeted high throughput parallel sequencing to detect the presence of specific mutations in RAS genes (KRAS, NRAS, HRAS, and BRAF). The test is performed on a 96-well plate and is suitable for use in a laboratory setting. The Praxis® Extended RAS Panel is a quantitative, colorimetric test using targeted high throughput parallel sequencing to detect the presence of specific mutations in RAS genes (KRAS, NRAS, HRAS, and BRAF). The test is performed on a 96-well plate and is suitable for use in a laboratory setting.	Количество: 1 шт. Объем: н/п Размеры: н/п	Описание: Инструкция по применению для набора реагентов Praxis Extended RAS Panel Бумага без покрытия белого цвета
Руководство справочное PN 15054440	 Praxis® Extended RAS Panel Reference Guide This reference guide provides detailed information about the Praxis® Extended RAS Panel, including its components, storage conditions, and handling instructions. It is intended for use by laboratory personnel who are responsible for the performance of the test. The Praxis® Extended RAS Panel is a quantitative, colorimetric test using targeted high throughput parallel sequencing to detect the presence of specific mutations in RAS genes (KRAS, NRAS, HRAS, and BRAF). The test is performed on a 96-well plate and is suitable for use in a laboratory setting.	Количество: 1 шт. Объем: н/п Размеры: н/п	Описание: Руководство справочное для Набора реагентов Praxis Extended RAS Panel Бумага без покрытия белого цвета

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IF exposed or concerned: Get medical advice/attention.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

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Praxis® Extended RAS Panel

1B

Library Normalization Diluent 1 x 4.6 ml

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Praxis® Extended RAS Panel

3/6

Stringent Wash Buffer	1 x 24 ml
Universal Wash Buffer	1 x 4.8 ml
PCR Clean-Up Beads	1 x 5 ml
Library Normalization Wash	2 x 4.8 ml
Library Beads	1 x 1.2 ml
MiSeqDx® Flow Cell	2 x 1

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DX# SW1

(Contains Formamide, Sodium Chloride)

May damage fertility.
May damage unborn child.
Causes serious eye irritation.

DANGER

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Use personal protective equipment as required. Wash thoroughly after handling. Wear protective gloves and eye/face protection.

If exposed or concerned: Get medical advice/attention.
If IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

www.illumina.com
1-800-809-4566

Label PN: 15041620.01

(gray shaded areas
represent varnish
over preprint)

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
Praxis® Extended RAS Panel

3B

PCR Clean-Up Beads	1 x 5 ml
Library Normalization Wash	2 x 4.8 ml
Library Beads	1 x 1.2 ml
MiSeqDx® Flow Cell	2 x 1

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DX LNWI

(Contains Citric Acid Monohydrate, Formamide, Sodium Chloride)

May damage fertility.
May damage unborn child.
Causes serious eye irritation.

DANGER

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Use personal protective equipment as required. Wash thoroughly after handling. Wear protective gloves and eye/face protection.

IF exposed or concerned: Get medical advice/attention.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

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


Praxis® Extended RAS Panel

5/6

Elution Buffer	1 x 4.8 ml
Library Storage Buffer	1 x 3.5 ml
TruSeq® Custom Amplicon Dx Filter Plate	2 x 1

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