Manufacturer:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

Manufacturing site:
Siemens Healthcare GmbH
Magnetic Resonance (MR)
Henkestr. 127
91052 Erlangen
Germany

User manual supplement for Russia The magnetic resonance system MAGNETOM Avanto with accessories

This document has been printed in Erlangen, Germany

A handwritten signature in blue ink, appearing to read "M. Velasco", written over a horizontal line.

Mayra Alejandra Velasco Alvarado
Regulatory Affairs
Siemens Healthcare GmbH, Magnetic Resonance (MR)

1. The name of the Medical Device: the accurate of the MD.

The magnetic-resonance system MAGNETOM Avanto with accessories.

2. Information about Manufacturer of the Medical Device:

2.1 Legal Manufacturer

Siemens Healthcare GmbH
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: siemens-healthineers.com

2.2 Manufacturing site

Siemens Healthcare GmbH - Magnetic Resonance (MR)
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: siemens-healthineers.com

2.3 The name and address of the local representative in Russia

Siemens Healthcare LLC
Address: Dubininskaya Str. 96, 115093 Moscow, Russia
Tel.: +7 495 737 12 52
Web: siemens-healthineers.com/ru/

3. Designation of the medical device: Purposes for use (Medical Device operation data):

3.1 Intended Use

The MAGNETOM Avanto system is indicated for use as a magnetic resonance diagnostic device (MRDD) that produces transverse, sagittal, coronal and oblique cross-sectional images, spectroscopic images and/or spectra, and that displays the internal structure and/or function of the head, body, or extremities. Other physical parameters derived from the images and/or spectra may also be produced. Depending on the region of interest, contrast agents may be used. These images and/or spectra and the physical parameters derived from the images and/or spectra when interpreted by a trained physician yield information that may assist in diagnosis.

The MAGNETOM Avanto system described above may also be used for imaging during interventional procedures when performed with MR compatible devices such as in-room display and MR Safe biopsy needles.

3.2 Indications

Please find following an overview of the indications based on available DOT Engine Imaging Protocols:

Dot Engine	Description
Brain Dot Engine	The Brain Dot Engine optimizes brain examinations with guided and automated workflows customized to your standards of care. The Brain Dot Engine supports the user in achieving reproducible image quality using automation tools and functionalities incorporated into the program.
Large Joint Dot Engine	The Large Joint Dot Engine optimizes image quality of knee, hip and shoulder scans by proposing the most appropriate protocols according to the examination strategy chosen for the specific patient. It ensures reproducible image quality and streamlines large joint

Dot Engine	Description
	examinations to the greatest extent. The Large Joint Dot Engine features AutoAlign and AutoCoverage for knee, hip and shoulder. The WARP and Advanced WARP techniques (including high bandwidth protocols, VAT and SEMAC) provide susceptibility artifact reduction functionality (e.g. from MR Conditional metal2) implants), and include optimized protocols for knee and hip examinations. High resolution 3D imaging programs together with user-configurable automatic Inline MPR (Multi Planar Reconstruction) calculations provide increased efficiency, reproducibility and ease of use.
Abdomen Dot Engine	The Abdomen Dot Engine offers standardized, efficient, and comprehensive workflows for the upper abdomen with excellent image quality. The workflow covers liver, biliary and pancreatic system and, if slightly adapted, kidneys as well. The workflow is prepared for easy reading and reporting together with syngo.via
Angio Dot Engine	The timing of contrast injection and scan is commonly stated as the most challenging part of an angiographic exam. The Angio Dot Engine guides the user through angiographic single or multi station examinations by providing visualization of arterial and venous timing windows using a test bolus technique. This information is fed back into the next planning steps so scan parameters can be adapted to the individual patient and patient's condition. Where needed, automatic voice commands support the communication with the patient.
Cardiac Dot Engine	Cardiac examinations used to be the most complex exams in MR. Now the Cardiac Dot Engine supports the user in many ways. Using anatomical landmarks, standard views of the heart, such as dedicated long axis and short-axis views are easily generated and can easily be reproduced using different scanning techniques. Scan parameters are adjusted to the patient's heart rate and automatic voice commands are given. All of this helps handle the complexity of CMR examinations with confidence and supports customized workflows that are easy to repeat.
Spine Dot Engine	The Spine Dot Engine delivers optimized cervical, thoracic and lumbar spine imaging for patients of all conditions and provides guided and automated workflows customized to your standards of care. The Spine Dot Engine supports the user in achieving reproducible image quality with increased ease-of-use and time efficient exams.
Breast Dot Engine	The Breast Dot Engine provides optimized protocols for tissue depiction, implant visualization, and breast biopsy. For ease of use, different examination strategies (FatSat, non-FatSat, feet-first/head-first positioning, InterVIEWS) are available for both high-channel coils and the 4ch BI coil.

3.3 Contraindications

In general, an MR examination is contraindicated for patients with electronic or electronically conductive implants or metals, especially those containing ferromagnetic foreign matter.

Typical contraindications for MR examinations are:

- Electronic implants: e.g. pacemakers, stimulators, insulin pumps
- Artificial heart valves, aneurysm clips
- Metal splinters in the eye (danger of retinal detachment)
- Artificial anus (anus praeter) with magnetic closure
- Transdermal drug patches with metallic backings
- Electrically conductive implants and prostheses
- Metallic spirals for contraception (IUDs = Intra Uterin Devices)
- Transdermal or other similar implants for body piercings as well as magnetic piercings

3.4 Possible side effects

Possible undesirable side effects for MR are dizziness, heating, claustrophobia and nerve stimulation.

3.5 Information about potential users of MAGNETOM Avanto

Typically physicians, trained radiological technicians, or technologists.

3.6 Mode of application (Medical Device application data)

1. Preparing an examination
 - Patient registration
 - Positioning the patient on the patient table
 - Positioning the local coils and cushions
 - Table positioning
2. Performing an examination
 - Measuring of the Localizer and parallel loading of previous reference images
 - Planning the examination volume
 - Scanning
 - Viewing and editing images during examinations
3. Basic syngo MR viewing, post-processing, and filming
 - Viewing
 - Mean Curve
 - Argus Viewer
 - Dynamic Analysis
 - 3D
 - Filming

3.7 Conditions for use

For ambient conditions in hospitals or MRI medical practices please refer as following:

Control room	
Temperature	15–30 °C (59–86 °F)
Relative air humidity	40 to 80%
Water content	11.5 g/kg dry air
Air pressure	1060–700 hPa ¹

Electronics room	
Temperature	15–30 °C (59–86 °F) (with integrated cabinet temperature control)
Relative air humidity	40 to 80%
Water content	11.5 g/kg dry air
Air pressure	1060–700 hPa ¹

Examination room	
Temperature	up to 22 °C
Relative air humidity	up to 60%
Water content	11.0 g/kg dry air
Air pressure	1060 – 7—hPa ¹

1. 700 hPa is equivalent to a height of about 3000 m above sea level

4. Classification of the Medical Device:

4.1 Type of the MD in accordance with the nomenclature for classification

- Europe (CE)

Kind of product:	Active medical device for diagnosis
Classification rules according to MDD 93/42/EEC, Annex IX:	Rule 10
Classification according to Annex IX:	Ila
Classification Criteria:	Active devices intended for diagnosis are in Class Ila if they are intended to supply energy which will be absorbed by the human body, except for devices used to illuminate the patient's body, in the visible spectrum
Classification in accordance with IEC 60601-1	Class I

- Russia

Classification in accordance with Russia Nomenclature Classification	135160
--	--------

4.2 Class of the potential risk of the medical device application

Potential risk in accordance to SW safety class "B": non-serious injury possible.

SW name and version (option 1)	syngo MR B19 Date of release: 28.06.2017
SW name and version (option 2)	syngo MR D13 Date of release: 11.07.2012
SW name and version (option 3)	syngo MR E11 Date of release: 02.09.2020

Please refer to chapter 5.4 for SW versions differences.

4.3 Type of contact with the human body

MR system device is not intended to be sterile.

4.4 Laser classification

The laser of the laser light localizer is classified as Class 2M according to IEC 60825-1.

4.5 IP Protection class

IP protection class according to IEC 60529 is IP X0.

4.6 Applied parts classification

Patient table	Type B
Local coils	Type B or BF
Body coil	Type B
ECG/Pulse module	Type BF

4.7 Operating mode

Continuous operation.


5. Technical Specifications of the Medical Device (the main technical parameters of the Medical Device; the MD's scheme (design), characteristics, properties, features):

5.1 The photo and flow chart of the device


Below you will find the front view and side view picture as well as the flow chart of the MAGNETOM Avanto.


Front view with fixed patient table - photo

Front view with fixed patient table –scheme in mm

Side view with fixed patient table - photo

Side view with fixed patient table –scheme in mm


Note: The cover of the system includes an official brand name of the system «MAGNETOM Avanto». Also the Cover of the system includes some of the following information:



- Tim [32 x 8].
- Tim [76 x 18].
- Tim [76 x 32].
- A Tim and Dot system

This additional information isn't a system's name and this information isn't indicated on the main system's labeling, User documentation etc. This additional information isn't part of the system's name and intends to be informative of the technology features of the system.



5.2 The list of components of the medical device


The magnetic-resonance system MAGNETOM Avanto with accessories			
##	Name of component	Description	Photograph
BASIC BLOCK:			
1.	MAGNETOM Avanto–System	Whole-body superconductive Zero Helium Boil-Off 1.5T magnet.	

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

3.	MAGNETOM Avanto – System	Integrated no tune transmit/receive Body Coil	
4.	Cover #Av T+D Cover #Avanto	Magnet cover	
5.	Patient Matrix Table #Av Remov. Matrix Table w.Trolley #Av 2nd Tabletop w/Trolley #Av	The patient table is mounted directly to the magnet assembly. The table comes with a comprehensive set of patient positioning cushions.	
6.	MAGNETOM Avanto – System	The electronics cabinet contains the electronics for Measurement sequence control generating the magnetic field gradients radio frequency generation Voltage supply for network Integrated cooling system The system shall support two cooling scenarios: passive and active cooling.	






7.	Separator #Av,Es	The SEP (Separation cabinet) has to be used if a central hospital chilled water supply is available or if a chiller of any brand/type is already available. In these cases, the primary water specifications must fulfill the requirements (e.g. 60kW heat dissipation; 90l/min flow; 6 to 12°C water temperature; ph value 6 to 8). Dimension: 1800 mm x 650 mm x 650 mm (height x width x depth) Weight: 400 kg Location: The position of the SEP has to be in the close vicinity of the ACC cabinet. Two 5m hoses (forward and return) to connect the SEP to the ACC, and two 1m hoses (forward and return) for connecting the SEP to the local water supply system are part of the delivery volume. Function: - Interface between the on-site water chiller (of any brand/type) or Interface to the central hospital chilled water supply. For the above-mentioned cases the SEP is mandatory! The cold head compressor is part of the SEP and is connected internal.	
----	------------------	--	--






8.	The alarm box	The alarm box has the following functions: - Displays alarm signals - Switches the MR system on and off - Magnet Stop/Magnet Quench The alarm box is installed near the syngo Acquisition Workplace.	
9.	Emergency switches	The MR system has different types of emergency switches. - Magnet Stop Shutting down the static magnetic field (quenching) - in case of accidents with attracted metal parts or in case of fire - Emergency Shutdown Electric power of the entire MR system is switched off, but magnet remains at field - E.g. in case of fire - Table Stop Motorized table movement is stopped in case of accidents or injury due to table movement.	
10.	The intercom	The intercom allows personnel and patients to communicate during the examination. In addition, some important operations like stopping the patient table can be managed from the intercom. Optionally, music or automatic voice outputs can be played in the examination room via the loudspeaker or the headphones.	 (1) Stop button (2) Talk button (3) Reset Patient Alert button (4) Listen button (5) Physio Signal button (6) Audio connection (7) Reset Stop button







11.	The syngo Acquisition Workplace	The syngo Acquisition Workplace (syngo Acq WP). It includes the host processor with the operating elements monitor, keyboard, and mouse.	
12.	PC Keyboard US English #Tim	The syngo Acquisition Workplace is equipped with an original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	
13.	Head Matrix Coil	Head Matrix Coil is used for applications like head examinations, MR Angiography, combined head/neck examinations, TMJ (temporo-mandibular joints)	




14.	Neck Matrix Coil	Neck Matrix coil is used for applications such as Cervical Spine, Neck, Larynx/Esophagus, MR Angiography, Mediastinum, combined head/neck examinations	
15.	Spine Matrix Coil	Spine Matrix coil is used for applications such as high resolution imaging of the whole spine, but also for various applications in combination with additional coils.	
16.	Body coil	The body coil with its 6 integrated pre-amplifiers ensures excellent signal-to-noise ratio. The SlideConnect plug allows for fast and easy patient preparation resulting in less table time.	
17.	Flex Coils	Light-weighted, very flexible, iPAT compatible, 4-element no-tune receiver coils which are made of soft and smooth material.	Flex Large 4 



			Flex Small 4 
18.	4-channel Flex Coil Interface	The Flex Coil Interface is required for connecting the Flex coils, the Loop coils and the Endorectal coil.	
19.	Connection cables	Connection cables	
20.	Measurement phantoms	Measurement phantoms are used for quality measurements.	Plastic bottle, 15ml  Plastic bottle, 1900ml  Plastic bottle, 5300ml  Spherical phantom D165  Spherical phantom D240 


			Hand/Wrist 16 phantom  Foot/Ankle 16 phantom  Special phantom for Endorectal 
21.	Phantom holder	Phantom holder	Phantom holder 17  Phantom holder 18  Phantom holder 23  Phantom holder 24 



			Phantom holder 25  Phantom holder 26  Phantom holder 27  Phantom holder 28  Phantom holder 29  Phantom holder 31 
--	--	--	--



			Phantom holder 32  Phantom holder 40  Phantom holder 45  Phantom holder 46  Phantom holder 47  Endorectal coil simulator 
--	--	--	--


			Prostate 2 phantom holder 
22.	The application package "Tim Application Suite": -Basic software license: basic license key for running the syngo software; -Basic MR applications license; -3D functionality; -Argus (cardiac) basic viewing; -DICOM worklist support; -PhoenixZIP; -Different standard measurement/examination and post-processing functionality.	The Tim Application Suite offers a complete range of clinically optimized sequences, protocols and workflow functionalities for all body regions.	Neuro Suite  Angio Suite  Cardiac Suite 




			Body Suite 
			Onco Suite 
			Ortho Suite 




			Breast Suite 
23.	Russian User Manual.	Instruction for use.	SIEMENS www.siemens.com/healthcare MAGNETOM Family Руководство пользователя – МР-система syngo MR E11 Answers for life SIEMENS www.siemens.com/healthcare MAGNETOM Family Руководство пользователя – МР-система syngo MR B19

			 MAGNETOM Avanto Инструкция по эксплуатации - МР-система серия MR D15 www.siemens.com/healthcare SIEMENS
	Принадлежности		
	syngo MR Workplace	The syngo MR Workplace allows for evaluation, documentation and post-processing of previously measured images while acquiring images at the syngo Acquisition Workplace.	
2.	Workstation-based Solution	syngo.via Workstation is a syngo.via-based entry level sales package for advanced visualization customers. It provides the same software functionality, as the conventional syngo.via product.	
3.	PC Keyboard US English	An original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	
4.	Foot switch	MR Foot Switch with two pneumatic buttons for start and stop of a preset MR sequence. The foot	



		switch is MR compatible and is positioned near the patient table on the examination room floor. The MR foot switch is especially suited for sequence control when MR examinations are carried out during operations.	
5.	PMU Wireless Physio Control #Ez	Physiological Measurement Unit (PMU) - Wireless Physio Control for wireless triggering, synchronizes the measurement with the physiological cycles of cardiac and/or respiratory motion.	 
6.	PMU Sensor Kit #Tim	Additional PERU, PPU and charging stand.	  


			 
7.	Respiratory cushion	Respiratory cushion with pressure hose	
8.	Respiratory belt	Respiratory belt	
9.	Application cushion	Application cushion	
10.	Charging station	The charging station is installed separately near the syngo Acquisition Workplace and is used for storing both units.	


11.	PMU Display #Av PMU Display Mobile #Av	The physiological signals are shown on the display	
12.	Control units	The control units are located at the right and the left side of the patient table at the front of the magnet cover.	Syngo MR B19  Syngo MRE11 
13.	Paper Roller Holder #Am	The Paper Roller Holder	 1 paper roll holder
14.	Headphones	The headphones are intended for communication with the patient and can be used in combination with ear plugs.	
15.	MagnaCoil	The MagnaCoil Headset can be used as an alternative headphone (for example, in head coils) to the standard Siemens headphones which are part of the products since years. Sometimes the standard headphones are not very practicable and therefore	



		Siemens wants to offer an alternative: the MagnaCoil Headset.	
16.	CP Head Array Coil	Dedicated coil for head examinations especially for utilization with a stereotactic frame for therapy planning and for pediatric imaging, e.g. of the spine. Especially for head exams with the use of a stereotactic frame.	
17.	32-channel Head Coil	The 32-channel Head Array Coil is an iPAT-compatible "no-tune" receiver coil for MAGNETOM Avanto	
18.	Tx/Rx CP Head Coil	Circularly polarized no-tune transmit/receive coil with an open patient-friendly design. The integrated transmit mode allows volume selective excitation. Integrated, extremely low-noise pre-amplifiers permit very high signal-to-noise ratio.	
19.	16-channel Body Coil, anterior #Tim 32-channel Body Coil #Tim	The Body Coil, anterior is the anterior part of the 32-channel Body Coil. It is an iPAT-compatible "no-tune" receiver coil for imaging of the heart /coronary arteries and imaging applications of the abdomen	16Ch Body Array Anterior MR Coil  16Ch Body Array Posterior MR Coil

			
20.	PA Matrix Coil #Tim	The PA Matrix Coil is also shipped with a set of positioning cushions for proper handling. No tuning of the fully iPAT-compatible Peripheral Angio 16 is required. Applications: - High-resolution angiography of both legs incl. pelvis (by additional use of the Body 13) with highest signal-to-noise ratio - Visualization of the iliac arteries and aorta in combination with Body 13 - Bilateral examinations of long bones of the legs.	
21.	4-channel Flex Coil, large #Tim	Light-weight, very flexible, iPAT compatible, 4-element no-tune receiver coil made of soft and smooth material. It is used for examinations of extremities (e.g. medium to large shoulder, hip or knee) and for abdominal examinations.	
22.	4-channel Flex Coil, small #Tim	Light-weight, very flexible, iPAT compatible, 4-element no-tune receiver coil made of soft and smooth material. It is used for examinations of the upper and lower extremities (e.g. small to medium shoulder, wrist, elbow or ankle).	

23.	Shoulder Array Coil #Av,S	This iPAT compatible coil for examinations of the left or right shoulder consists of a base plate and two receive array coil attachments available in different sizes, these will be attached and can be relocated on the basis plate.	
24.	8-channel Knee Coil #Tim	Knee Coil is an iPAT compatible no tune receive coil for knee joint examinations.	
25.	Tx/Rx 15-channel Knee Coil #Tim	Coil for joint examinations in the area of the lower extremities.	
26.	8-channel Wrist Coil #Tim	The wrist coil is an iPAT-compatible "no-tune" receiver coil for the examination of the wrist and fingers.	

27.	8-channel Foot/ Ankle Coil #Av,S	The foot-ankle coil is an iPAT-compatible "no-tune" receiver coil for the examination of the foot and the ankle joint.	
28.	CP Extremity Coil #Tim	High-performance extremity coil, optimized for knee imaging.	
29.	Breast Matrix Coil #Tim CP Breast Array Breast Array	The Breast coil delivers brilliant image quality for high-resolution 2D and 3D MR breast imaging. Compression devices adjustable in increments enable stabilization of the breast during the examination. In addition, anatomically shaped cushions and pads are optionally supplied for optimal patient positioning.	
30.	2/4/8-ch Sentinelle Breast 2/10/16-ch Sentinelle BreastCoil	The 16-channel imaging configuration of the Sentinelle Breast Coil consists of two lateral 4-channel coil elements and an 8-channel coil middle element. The 16-channel Sentinelle Breast Coil delivers brilliant image quality for high-resolution 2D and 3D MR breast imaging. Techniques for reducing scan times, such as	2-/4-/8-Channel Sentinelle Breast Coil  2-/10-/16-Channel Sentinelle Breast Coil 



		parallel imaging, can be used very well. The biopsy configuration consists of two lateral 1-channel coil elements and an 8-channel coil middle element. For the unilateral biopsy setup a contralateral support will be used. The Sentinelle Breast Coil supports the Grid biopsy method.	
31.	Positioning frame f. BreastCoil#T+D	Open design positioning frame for inserting the 4-channel BI Breast Coil. The positioning frame comes supplied with: - a head cushion - fixation rail to proper position and fix the coil positioning frame on the patient table - abdominal support cushions - Breast stabilizing cushions - positioning guidance	
32.	16-channel AI Breast Coil	The 16-channel AI Breast Coil consists of a 16-element design with 16 integrated preamplifiers in which the coil elements are arranged on each side in an array of 6 plus an axilla element and a surrounding loop element. The 16-channel Breast Coil delivers brilliant image quality for high-resolution 2D and 3D MR breast imaging. The 16-channel AI Breast Coil measures approx. 71 mm x 47 mm x 20 mm (LxWxH) and weighs approx. 10 kg with positioning frame.	
33.	4-channel BI Breast Coil #Tim 7 Ch Breast	The coil consists of a 2-channel insertion plate and a 2-channel support coil on each side.	Insertion Plate (left) Patient Pad (right)




			
34.	BI/AI Coil package + Frame	The package consisting of the 4-channel BI Breast Coil, Biopsy Training Starter Kit and the 16-channel AI Breast Coil is an attractive option if a breast biopsy coil with a large biopsy access and an outstanding imaging coil offering 16 channels is needed. The package includes a frame in which the lower part and the upper part of the 4-channel BI Breast Coil can be placed. The same frame is also suitable for use with the 16-channel AI Breast Coil.	
35.	Special Purpose Coil #Tim	The 4-channel special purpose coil is a no-tune receiver coil designed for examination with small Field-of-Views.	
36.	Loop Coil, large Loop Flex Coil, large #4-ch	No-tune loop receive coil for examinations of the upper or lower extremities (e.g. shoulder).	Loop Coil Large  Flex Loop Large 4Ch


			
37.	Loop Coil, medium	No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin).	
38.	Loop Coil, small	No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin).	
39.	Double Loop Array Coil	The Double Loop Array Coil consists of two flexible, anatomically adaptable 7 cm ring coils for simultaneous examination of both TMJs, the eyes or the wrists with optimized, excellent image quality.	
40.	Biopsy Starter Kit	The breast biopsy kit comes as a case containing a fixation unit plus accessories for both the Post-Pillar and the Grid biopsy method. This is a training kit, not used with patients.	Case 




			Post/Pillar Grid
41.	Endorectal Coil Kit	Kit including interface device for connecting the prostate receive coil probe. The interface device connects to the Tim Coil Interface.	
42.	Prostate 2	The Prostate 2 is used for examinations of the prostate, colon and the surrounding tissue in the pelvis. Depending on the requirements of the examination, the coil may be used in the Head First or Feet First orientation.	
43.	Prostate Endo.Coil	The Medrad disposable Prostate coil, through its high-resolution imaging and small field of view, provides improved ability to visualize internal architecture of the prostate and periprostatic structures,	





		including prostate capsule and vascular bundles. The radiological images may be used in diagnosing diseases of the prostate. The inflated balloon makes a positioning of the Prostate coil in close proximity to the wall adjacent to the prostate gland possible.	
44.	Colorect. Endo.Coil	The Medrad disposable Colorectal coil, through its high-resolution imaging and small field of view, provides images of the colorectal anatomy, which may be used to aid in diagnosing colorectal diseases.	
45.	Cervix Endo.Coil	The Medrad disposable Cervix coil, through its high-resolution imaging and small field of view, provides images of cervix, uterine, and endometrial anatomies, which may be used to aid in diagnosing diseases in these anatomical regions.	
46.	Pediatric coil	A head and neck coil with excellent SNR for small FOV, high resolution, fast imaging of children up to 18 months old. The open design allows easy access for intubation and observation of the patient. The infant cradle can be used to help transport and or position very small patients and is recommended for children no older than 5 months of age.	



47.	31P/1H Heart-/ Liver Coil #S,Av	Dual-resonant 1H/31P surface coil for heart and liver spectroscopy.	
48.	4-channel Flex Coil Interface #Tim	Interface with integrated preamplifiers for the connection of the following flexible coils to the system: - Flex Large 4 - Flex Small 4 - Loop Coil, large - Loop Coil, medium - Loop Coil, small - Endorectal Coil	
49.	Tim Coil Interface 1.5T Adapt.Sy/So/Es Should.Arr.Coil#Av,S 4+1 Adapter #S,Av,Es	Coil adapter plug for up to 8 receive and 1 transmit channels, in order to connect existing dedicated knee and breast coils (Tx/Rx 15-channel Knee Coil, CP Extremity Coil, 4-channel BI Breast Coil, 16-channel AI Breast Coil, (2/4)/8-channel Sentinelle BreastCoil and (2/10)/16-channel Sentinelle BreastCoil) with all MAGNETOM 1.5T Systems using Tim 4G-technology.	
50.	Coil Storage Cart #T+D,Ez	Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories.	





51.	Accessoires cabinet	Two designer cabinets matching the magnet design, thus enhancing the patient-friendly impression. The cabinets are done in the colors of the magnet and serve the purpose of accommodating the coils and examination accessories.	
52.	Cushions and positioning aids	Cushions and positioning aids - Table pad 1 - Table pad 2 - Table pad 3 - Table pad 4 - Table pad 5 - Head insert - 1 paper roll holder - 4 belt straps - 2 different-sized sandbags, 2 each	 Table pad 1  Table pad 2  Table pad 3  Table pad 4  Head-heel insert





			 1 paper roll holder  4 belt straps  2 different-sized sandbags, 2 each
53.	Set of positioning aids	Set of positioning aids - 1 knee wedge - 1 positioning roll - 1 back-of-the head support - 2 large triangles - 2 small triangles - 1 wedge - 1 lordosis cushion - 2 flat cushions - 2 flat wedges - 2 knee cushions - 2 different-sized multipurpose cushions, 3 pieces each - Gap cushion set small - Gap cushion set large	 1 knee wedge  1 positioning roll  1 back-of-the head support




			 2 large triangles  2 small triangles  1 wedge  1 lordosis cushion  2 flat wedges  2 knee cushions
--	--	--	---




			 2 different-sized multipurpose cushions, 3 pieces each  Gap cushion set small  Gap cushion set large
54.	Head Matrix Coil	Head Matrix Coil accessories - Relaxation cushion - Head support - Head stabilization wedge - Mirror carrier	 




			 
55.	CP Head Array	CP Head Array accessories - Thick cushion - Thin cushion - Head support	 





			
56.	8ch Head Array	8ch Head Array accessories - 2 blue positioning aids	
57.	TxRx CP Head	TxRx CP Head accessories - Head cushion - 2 head wedges (thin size) - Head-Spine table pad - 2 head wedges (thic size)	 Head cushion  2 head wedges (thin size)



			 Head-Spine table pad  2 head wedges (thick size)
58.	TxRx Head	TxRx Head accessories - Head cushion - 2 head wedges	 
59.	32 Ch Head	32 Ch Head accessories - Mirror carrier with mirror - Cushion 32 - Head alignment wedges	





			 
60.	Double Loop Array	Double Loop Array accessories - Hand rest - Thick cushion	 




61.	Neck Matrix	Neck Matrix accessories - Neck cushion - 2 positioning aids for the neck	 
62.	Spine Matrix	Spine Matrix accessories - Cushion support	 Cushion support
63.	Body Matrix	Body Matrix accessories - iPAT cushion - Elastic connection clamps - 3 cable holders - 2 straps - 6 belts - 4 long straps	



			   
--	--	--	--


			
64.	16 Ch Body Array A/P	16 Ch Body Array A/P accessories - Base positioner - 2 Comfort pads 150 - 2 Comfort pads 430 - Patient support pad - Velcro straps	  





			 
65.	Heart/Liver	Heart/Liver accessories - Positioning pad - Heart/Liver coil cushion - Additional Heart/Liver cushion - 2 straps - 4 belt straps - Straps	 


			   
--	--	--	--





66.	Endorectal Coil	Endorectal Coil accessories - Endo cushion - Endorectal Adapter (1.5T coil)	 Endo cushion  Endorectal Adapter 1.5T
67.	Prostate 2	Prostate 2 accessories - Patient support cushion - Back pad - Foot rest pad - Body Matrix cushion - Plug cap - Pressure testing cap - Cable extension - Base plate - Pressure testing unit	 Patient support pad  Back pad  Foot rest pad  Body 18 cushion





			 Plug cap  Pressure testing cap  Cable extension  Base plate  Pressure testing unit
68.	Breast Matrix	Breast Matrix accessories - 2 large compression attachments - 2 small compression attachments - Positioning aids - Soft pad - Head support - Wedge cushion (abdominal support) - Head cushion - Support for cushion - Head pad	









			   
--	--	--	--




69.	CP Breast Array	CP Breast Array accessories - 2 large compression attachments - 2 small compression attachments - Positioning aids - Pad - Head support - Wedge cushion (abdominal support)	
-----	-----------------	--	---




			 
70.	Breast Array	Breast Array accessories - Pad - Head support - Wedge cushion (abdominal support)	 




			
71.	7 Ch Breast	7 Ch Breast accessories - Pad - Head support	 
72.	Shoulder Array	Shoulder Array accessories - Special head support - Pad	




			
73.	8 Ch Wrist	8 Ch Wrist accessories - Base pad - Elbow pad small - Finger pad - Palm pad - Upper pad - Wrist pad - Elbow pad thin - Elbow pad	  




			   
--	--	--	--



			
74.	PA Matrix	PA Matrix accessories - Belt pad (1.5T coil) - Additional cushion for large patients (1.5T coil) - 2 knee cushions - 4 flat wedges - 2 flat cushions	  



			 
75.	CP PA Array	CP PA Array accessories - Belt pad - 2 knee cushions - 4 flat wedges - 2 flat cushions	 



			 
76.	High Resolution Knee TxRx 15 Ch Knee	High Resolution Knee accessories - HRK Coil Pad, 0,25 - HRK Coil Pad, 0,50 - HRK Ankle Pad - HRK Ramp Pad	 


			 
77.	8 Ch Knee	8 Ch Knee accessories - HRK Coil Pad, 0,25 - HRK Coil Pad, 0,50 - HRK Ankle Pad - HRK Ramp Pad	- HRK Coil Pad, 0,25;  - HRK Coil Pad, 0,50;  - прокладка для голеностопа HRK;

			 - пандусная подушка HRK. 
78.	CP Extremity	CP Extremity accessories - Knee joint cushion - Leg support - Large ankle cushion - Small ankle cushion - Foot cushion - Cushion 4 - Foot support	 Knee joint cushion  Leg support  Large ankle cushion


			 Small ankle cushion  Foot cushion  Air cushion  Foot support
79.	8 Ch Foot Ankle	8 Ch Foot Ankle accessories - Knee cushion - Heel cushion - Foot cushion	 Knee cushion  Heel cushion




			 Foot cushion
80.	Flex Small/Large 4	Flex Small/Large 4 accessories: - Hand wedge - 2 soft straps - 3 cable holders - 4 belts - iPat cushion, small - iPat cushion, large	 Hand wedge  2 soft straps  4 belts  iPat cushion, large  iPat cushion, small


81.	Comfort Kit Vacuum cushion	The kit consists of an integrated vacuum pump and 3 anatomically shaped (spine, head, small joints) cushions of different sizes: - Head cushion - Spine cushion - Small cushion.	 Head cushion  Spine cushion  Small cushion
82.	Optional Dot Engine Workflows	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
83.	Optional Application Packages	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	

			
84.	Optional Postprocessing Packages	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
85.	Additional optional single Applications	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
86.	Argus Flow Argus Dynamic Signal	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	

87.	syngo Expert-i syngo Expert-i MRWP #T+D	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
88.	Syngo Security package	Software, optional application package. Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. This package supports customers in their achieving compliance with the HIPAA "Privacy" rule.	
89.	RF-Cabin #1.5T RF-Cabin #Ae,Av,Am,Se	Steel-full-room RF cabin, in a module method of building. Steel-full-room RF cabin, in a module method of building.	
90.	RF-Traverse Beams #1.5T RF-Traverse Beams #Ae,Av,Am,Se	4 traverse beams for the RF cabin. 5 traverse beams for the RF cabin.	

91.	RF-Interior #1.5T RF-Interior #Ae,Av,Am,Se RF-Interior, LED #1.5T RF-Interior, LED #Ae,Av,Am,Se	Standard interior décor with Halogen lighting for the RF cabin. Standard interior décor with Halogen lighting for the RF cabin. Standard interior décor with LED lighting for the RF cabin. Standard interior décor with LED lighting for the RF cabin.	
92.	RF Sound Proofing #Av,Es,Ez,S RF Sound Proofing #Ae,Av,Am,Se	Basic sound protection system.	
93.	Filter Cabinet, big Filter Cupboard	Cabinet for filter plate inside of the RF cabin. Wall-mounted cabinet custom-made. Cupboard for filter plate inside of the RF cabin. Wall-mounted cabinet custom-made.	
94.	Filter plate for RF-Cabin	Filterplate in low leakage circuitry for shielded rooms (MRI) 6 separate lines, 250 V, 50/60 Hz, 250Vdc, 32A	
95.	Adapter, RF Cabin #Tim,Ez	2 adapter plates for adjusting the larger filter plate opening of an	Filter plate with cabling:

		existing RF-room enclosure.	 
96.	Threshold-free door	Soundproof threshold-free door for the RF-cabins that can be opened and closed without effort by pneumatic mechanism.	
97.	Emergency hatch	With decorative plate on both sides, white, to be opened from the outside, approx. 600 x 600 mm.	



98.	Observation/ Day-light window	Observation / daylight window that can be opened, installation sideways.	
99.	Patient Supervision TV #Av,TATS Patient TV wall support	This package contains a special video camera for monitoring the patient during an MR examination, conveniently mounted on the wall of the examination room. The supervision solution is customizable and designed to address different site specific requirements. Up to 4 cameras can be optionally connected for patient supervision in the examination or waiting room. Additional video camera to allow for comprehensive monitoring of the patient's face independent of patient position during an MR examination. Wall mount for the patient monitor.	


100.	Video Distributor for monitor #T+D	Video distributor for connection of up to two additional color LCD monitors as well as an In-Room Acq WP syngo (in the examination room parallel to the syngo Acquisition Workplace operator monitor).	
101.	LCD Color Remote Monitor #1.5T	High-resolution LCD color flat screen monitor, with integrated gamma correction for optimum display of MR grayscale images and automatic backlight control for long-term brightness stability.	
102.	Color monitor	Color monitor	
103.	Trigger Converter	An optical trigger signal is available to trigger external stimulation devices in fMRI experiments.	Trigger Converter   
104.	MR Workplace Table	Ergonomically designed table for : - Color monitor - Keyboard - Mouse - Patient communication unit - Patient supervision display	 (1) MR Workplace Table (2) MR Workplace Container, 50 cm (optional)

105.	MR Workplace Container	Extra case for the syngo host computer with sliding front door to allow change of storage media	 (1) MR Workplace Table (2) MR Workplace Container. 50 cm (optional)
106.	MRI-Stretcher, non-magnetic MRI-Wheelchair, non-magnetic	MRI-Stretcher, slightly magnetic, chrome-nickel compounded for patient transport inside the clinic building. Continuously adjustable head piece with spring brake, synthetic leather cushion, movable barrier on either side, stretcher movable on four fixable wheels diameter 100 mm.	
107.	IV Stand, non-magnetic	IV-Rod, fix mounted, for AGA non magnetic stretcher. - Diameter: 180mm - Length:1720mm - With four hooks	
108.	Glider Positioning Kit	Positioning System for patient transport between bed or wheelchair and MR-System. These sliding boards consist of a special material with very low coefficients of friction. They can be folded, take up little space and can be transported easily. They are firm but flexible and adapt to different shapes of the body.	

109.	Instrument table, non-magnetic	Mobile 4-leg frame construction made of steel. Two shelves: surface ground to reduce glare, resistant to disinfectants, removable. 4 plastic swivel castors, two of them with brakes, wear resistant, non-discoloring when in contact with other material, electrically conductive.	
110.	Upgrade Magnet Cooling System	Exchange of the magnet cooling components (helium compressor and cold head) by models of the latest generation.	
111.	UPS system UPS Battery module	UPS system Liebert GXT4 3000RT230E for safeguarding computers. Including Power Cable of 9 m for connecting the UPS. UPS battery module Liebert GXT4 72VBATTE for safeguarding computers.	

112.	Distribution board	Distribution board	
113.	Metal detector	Metal detector	
114.	A fire extinguisher	A fire extinguisher	

115.	Seismic Kit	It comprises: - 3/4" tufnol spacers which have hole offsets from zero to 4.2 mm in 0.7 mm increments and washers - Radial gap checking tool - Vertical gap setting tool - 1" tufnol spacer kit - only to be used with 1" floor anchors - Loctite 243 (10 ml bottle) - 4 x M20 washers	
116.	ECG-electrods	ECG-electrods	
117.	Service Syphon Kit	Service Syphon Kit	

118.	HE-SYPHON DEWAR END	HE-SYPHON DEWAR END	
119.	SYPHON EXT.LEG	SYPHON EXT.LEG	

120.	Floor plate for table Trolley Heightening Kit #Av,Es,SATS	Base Plate for the fixed table	
121.	Adhes. Plastic-Mastic complete	Adhesive Plastic-Mastic complete kit consisting of adhesive component A, adhesive component B, gloves and spatle.	

5.3 Technical requirements

Technical Requirements based on the key parameters and values as presented in the chapter below.

5.4 Key parameters & values provided by manufacturer

- syngo MR B19 technical data:

Magnet parameters	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Number of field generating coils	7
Magnet length	150 cm
Magnet bore diameter ¹	90 cm
Inner tunnel length	88 cm
Length of flared opening	54.5 cm
Inner diameter ²	60 cm
System length	160 cm
Type of installation	Fixed
Decay characteristics from full field to 20 mT	approx. 20 seconds

1. without shim coils, gradient coil, RF body coil

2. incl. shim coils, gradient coil, RF body coil

Homogeneity (based on highly accurate 24 plane plot)		
DSV	Guaranteed	Typical
50 cm	<1.5 ppm	0.8 ppm
45 cm	<1 ppm	0.4 ppm
40 cm	<0.75 ppm	0.2 ppm
30 cm	<0.5 ppm	0.1 ppm
20 cm	<0.25 ppm	0.04 ppm
10 cm	<0.05 ppm	0.01 ppm

DSV – Diameter spherical volume (x, y and z direction)

Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Standard active shim with 3 linear channels

3D Shim	
• Patient-specific automated shim • For optimization of the magnetic field homogeneity	
Time to shim	approx. 20 s

Active Shielding (AS)		
Fringe field (axial x radial)	0.5 mT ¹	4.0 m x 2.5 m
	0.1 mT	6.0 m x 3.5 m

1. pacemaker safety limit

External Interference Shield (E.I.S.)
Patented shielding system integrated into the magnet
Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects)

Magnet cooling system	
Refill interval (typical)	10 years ¹
Boil-off rate (typical)	0.0 l/year ¹
Helium capacity	approx. 1600 liters
Minimum helium level	40%
Cryostat	Stainless steel

1. For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

Patient handling:

The patient table is available in two different configurations:

- Matrix configuration (standard)
- Whole-body configuration (Tim Whole Body Suite option) The respective configuration is set by Siemens Service during MR system installation.

Patient table	
Max. patient weight (including vertical movement)	250 kg
Max. scan range	154 cm optional 205 cm ¹
Vertical table movement	
Range	47–89 cm
Max. speed	10 cm/s
Horizontal table movement	
Max. range	263.5 cm
Max. speed	20 cm/s
Horizontal accuracy for repositioning from one direction	±0.8 mm

1. with Tim Whole Body Suite option

Patient comfort	
Inner diameter of magnet bore ¹	60 cm
In-bore lighting	can be set at 3 different levels
In-bore ventilation	can be set at 3 different levels
In-bore intercom	including loudspeaker, microphone and earphones

1. incl. shim coils, gradient coil, RF body coil

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2010 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise	
Patients require hearing protection with an SNR ¹ = 11 dB or more	

1. SNR = Single Number Rating

Noise for personnel in the examination room	
Noise measured	75.0 dB (A)

Gradient system

General features	
Actively shielded (AS), super-compact whole-body gradient coil system	
Extremely low eddy currents	
Water-cooled coil and amplifier for maximum performance	
Gradient coil length	122 cm

Patient-specific shimming	
Active shim via 3 linear channels with 20 coils	
Prepared for 5 additional non-linear 2nd order channels with 32 coils for advanced high-order shim (optional)	

Gradient amplifier		
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses		
Q-engine	Max. output voltage ¹	1250 V
	Max. output current ¹	460 A
SQ-engine	Max. output voltage ¹	2000 V
	Max. output current ¹	625 A

1. Values for each of the 3 gradient axes

TIm [32×8] Q-engine	
Performance for each axis	
Max. amplitude	33 mT/m
Min. rise time	264 µs from 0 to 33 mT/m
Max. slew rate	125 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s

TIm [76×18] Q-engine	
Performance for each axis	
Max. amplitude	33 mT/m
Min. rise time	264 µs from 0 to 33 mT/m
Max. slew rate	125 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s

T1m [76×18] SQ-engine	
Performance for each axis	
Max. amplitude	45 mT/m ¹
Min. rise time	200 µs from 0 to 40 mT/m
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	72 mT/m
Max. eff. slew rate	346 T/m/s

1. Value is valid for longitudinal direction. Max. amplitude for horizontal and vertical direction: 40 mT/m

T1m [76×32] SQ-engine	
Performance for each axis	
Max. amplitude	45 mT/m ¹
Min. rise time	200 µs from 0 to 40 mT/m
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	72 mT/m
Max. eff. slew rate	346 T/m/s

1. Value is valid for longitudinal direction. Max. amplitude for horizontal and vertical direction: 40 mT/m

RF system

General	
• Solid-state amplifier, maintenance-free • Digital quadrature demodulation • Digital filtering • Receiver with high dynamic range without adjustments	
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$
Frequency control	32 bits (0.015 Hz)
Phase control	16 bits (0.006 degrees)

Transmitter path	
Transmit amplitude	16 bit control 50 ns resolution
Gain stability (5 min)	<0.1 dB (low signal path) <0.4 dB (power amplifier)

Transmit amplifier	
Extremely compact water-cooled	
Bandwidth	500 kHz
Peak power	15 kW

Receive path	
Bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bits
Number of independent receiver channels	8/18/32
Sampling rate (Sampling resolution)	10 MHz (100 ns)
Preamplifier noise figure	<0.5 dB
Total gain	36 dB/51 dB with automatic control
Dynamic range (referred to 1 Hz resolution bandwidth)	162 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

Multi-Nuclear Option¹				
• Multi-Nuclear Imaging and Spectroscopy for the following nuclei: 3He, 7Li, 13C, 19F, 23Na, 31P, 129Xe • Decoupling nuclei: 13C, 23Na, 31P • Hyperpolarization nuclei: 3He, 13C, 129Xe				
Frequency band	3He	7Li	13C	19F
Center frequency [MHz]	48.4498	24.7172	15.9933	59.8336
Usable transmitter band-width [kHz]	±50	±50	±50	±50
Usable receiver band-width [kHz]	±200	±105	±70	±50
	23Na	31P	129Xe	
Center frequency [MHz]	16.8233	25.7455	17.5941	
Usable transmitter band-width [kHz]	±50	±50	±50	
Usable receiver band-width [kHz]	±70	±105	±35	

1. Optional

Computer system

Host computer (syngo Acquisition Workplace)	
Processor	1×Intel® Xeon® Quad-Core
Clock rate	≥3.5 GHz
Main memory (RAM)	8 GB
1 st hard disk (system SW)	300 GB
2 nd hard disk (data base)	300 GB
3 rd hard disk (images)	300 GB (approx. 110000 images, matrix 256 ² and 512 ² uncompressed)
CD-R writer	approx. 4000 images 256 ² ; DICOM Standard, ISO 9660
DVD-R writer	approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
Media drives	CD/DVD drive
Color LCD monitor ¹	Screen matrix 1280×1024 pixels

1. A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

Image Processor	
Processor	≥1 Intel® Xeon® Quad-Core
Clock rate	≥3.3 GHz
Main memory (RAM)	≥16 GB
1 hard disk (system SW)	240 GB SSD
4 hard disks (raw data)	240 GB SSD
Reconstruction speed	≥8511 recons per second (256 ² FFT, full FoV) ≥40816 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 8 data sets
Media drives	CD/DVD drive

syngo MR Workplace (optional)	
Processor	1×Intel® Xeon® Quad-Core
Clock rate	≥3.5 GHz
Main memory (RAM)	8 GB
1 st hard disk (system SW)	300 GB
2 nd hard disk (images)	300 GB
CD-R writer	approx. 4000 images 256 ² ; DICOM Standard, ISO 9660
DVD-R writer	approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
Media drives	CD/DVD drive
Color LCD monitor ¹	Screen matrix 1280×1024 pixels

In-Room syngo Acquisition Workplace (optional)	
Color LCD monitor ¹	
Screen size (diagonal)	18"
Screen matrix	1280×1024 pixels

1. A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

Euro Connector (System IEC320)	
Load rating	100–240 Vac
	50–60 Hz
	max. 10 A

Network/data coupling

The MR system provides the level of safety according IEC 60950-1 outside patient environment. All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

Cooling unit

Water consumption	Water cooled cabinet	>90 l/min ¹
Heat dissipation to water	Water cooled cabinet	
	Maximum	SQ-engine ≤60 kW Q-engine ≤45 kW
	Ready for measurement ¹	≤11.5 kW
	Typical measurement ¹	SQ-engine ≤16 kW
		Q-engine ≤15.5 kW

1. Water temperature: 6–12 °C (43–54 °F)

Line power supply

Variable	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50 Hz, 60 Hz	±1 Hz
Connection voltage ¹	Q-engine: 70 kVA SQ-engine: 110 kVA	
On-site protection	at least 160 A	

1. Long-term average in full operation

For all products, line power has to be supplied via an on-site system contactor or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710.¹

1. In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plate

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT – TIME	kVA: 105
LONG – TIME	kVA: 25

Power consumption

System off ¹	≤5.0 kW
Stand-by ¹	≤5.7 kW
Ready for measurement ¹	≤11.5 kW
Typical measurement ¹	Q-engine: ≤15.5 kW
	SQ-engine: ≤16.0 kW

1. All data incl. cold head compressor, without chilling

Dimensions

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipation [kW]
Examination Room					
Magnet 1.5 Tesla AS (incl. Helium)	207	150	210	4050±80	
Magnet in operation, incl. gradient coil, body coil, patient table, and covers	229.5	405	230	5850	3
Patient table	54 (top)	243	47–89	350	
Required min. room height clearance			235		
Min. transport dimensions	230/207 ¹	305/170 ¹	210 ¹		
Control Room					
syngo Acquisition Workplace (table+monitor)	120	80	117 (72+45)		2
Host computer	22	46	47		
syngo MR Workplace (table+monitor)	120	80	117 (72+45)		
Equipment Room					
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	160	65	189 ²	1500	5
Cooling system	65	65	189 ²	500	

1. without covers, electronics, table and table bracket

2. without attachments

- syngo MR D13 technical data

Magnet system

Magnet parameters	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Number of field generating coils	7
Magnet length	150 cm
Magnet bore diameter ¹	90 cm
Inner tunnel length	88 cm
Length of flared opening	54.5 cm
Inner diameter ²	60 cm
System length	160 cm

1. without shim coils, gradient coil, RF body coil

2. incl. shim coils, gradient coil, RF body coil

Homogeneity (based on highly accurate 24 plane plot)		
DSV	Guaranteed	Typical
50 cm	<1.5 ppm	0.8 ppm
45 cm	<1 ppm	0.4 ppm
40 cm	<0.75 ppm	0.2 ppm
30 cm	<0.5 ppm	0.1 ppm
20 cm	<0.25 ppm	0.04 ppm
10 cm	<0.05 ppm	0.01 ppm

DSV = Diameter spherical volume (x, y and z direction)

Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Standard active shim with 3 linear channels

3D Shim	
• Patient-specific automated shim • For optimization of the magnetic field homogeneity	
Time to shim	approx. 20 s

Shielding		
Active Shielding (AS)	5 th generation active shielding (AS) technology with counter coils	
Fringe field (axial×radial)	0.5 mT ¹	4.0 m×2.5 m
	0.1 mT	5.5 m×3.1 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

1. pacemaker safety limit

Magnet cooling system	
Refill interval (typical)	Not applicable ¹
Boil-off rate (typical)	0.0 l/year ¹
Helium capacity	approx. 1600 liters
Minimum helium level	40%
Cryostat	Stainless steel

1. For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

Patient handling

Tim table		
Max. patient weight (including vertical movement)		250 kg
Max. scan range		154 cm optional 205 cm ¹
Vertical table movement	Range	47–89 cm
	Max. speed	10 cm/s
Longitudinal table movement	Max. range	263.5 cm
	Max. speed	20 cm/s
	Positioning accuracy	±0.8 mm

1. with Tim Whole Body Suite option

Patient comfort	
Inner diameter of magnet bore ¹	60 cm
In-bore lighting	can be set at 3 different levels
In-bore ventilation	can be set at 3 different levels
In-bore intercom	including loudspeaker, microphone and earphones

1. incl. shim coils, gradient coil, RF body coil

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2006 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise	
Patients require hearing protection with an SNR ¹ = 11 dB or more	

1. SNR – Single Number Rating

Noise for personnel in the examination room	
Noise measured	75.0 dB (A)

Gradient system

General: actively shielded (AS) whole-body gradient coil system

Gradient amplifier		
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses		
Q-engine	Max. output voltage ¹	1250 V
	Max. output current ¹	460 A
SQ-engine	Max. output voltage ¹	2000 V
	Max. output current ¹	625 A

1. Values for each of the 3 gradient axes

Tim [76×18] Q-engine	
Performance for each axis	
Max. amplitude	33 mT/m
Min. rise time	264 µs from 0 to 33 mT/m
Max. slew rate	125 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s

Tim [76×18] SQ-engine	
Performance for each axis	
Max. amplitude	45 mT/m ¹
Min. rise time	200 µs from 0 to 45 mT/m
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	72 mT/m
Max. eff. slew rate	346 T/m/s

1. Value is valid for longitudinal direction. Max. amplitude for horizontal and vertical direction: 40 mT/m

Tim [76×32] SQ-engine	
Performance for each axis	
Max. amplitude	45 mT/m ¹
Min. rise time	200 µs from 0 to 45 mT/m
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	72 mT/m
Max. eff. slew rate	346 T/m/s

1. Value is valid for longitudinal direction. Max. amplitude for horizontal and vertical direction: 40 mT/m

RF system

Digital signal processing system	
• Solid-state amplifier, maintenance-free • Digital quadrature demodulation • Digital filtering • Receiver with high dynamic range without adjustments	

RF transmit technology		
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.015 Hz)	
Phase control	16 bits (0.006 degrees)	
Transmitter path	Transmit amplitude	16 bit control 50 ns resolution
	Gain stability (5 min)	<0.1 dB (low signal path) <0.4 dB (power amplifier)
Transmitter amplifier	Extremely compact water-cooled	
	Transmit amplifier band-width	500 kHz
	Peak power	15 kW

Receive path	
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bits
Number of independent receiver channels	18/32
Sampling rate (Sampling resolution)	10 MHz (100 ns)
Preamplifier noise figure	<0.5 dB
Total gain	36 dB/51 dB with automatic control
Dynamic range (referred to 1 Hz resolution bandwidth)	162 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

Computer system

syngo Acquisition Workplace		
Host computer	Processor	Intel Xeon ≥ W3520 Quad-Core
	Clock rate	≥2.66 GHz
	Main memory (RAM)	≥6 GB
	1 st hard disk (system SW)	300 GB SAS
	2 nd hard disk (data base)	300 GB SAS
	3 rd hard disk (images)	300 GB SAS
	CD-R writer	Approx. 4000 images 256 ² ; DICOM Standard, ISO 9660
	DVD-R writer	Approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
	Media drives	CD/DVD drive
Color LCD monitor	Screen matrix	1280×1024 pixels
Measurement and reconstruction system	Processor	Intel ≥ E3-1225 Quad-Core
	Clock rate	≥3.1 GHz
	Main memory (RAM)	16 GB for Tim [76×18] 32 GB for Tim [76×32]
	Hard disk for raw data	≥500 GB
	Hard disk for system software	≥500 GB
	Reconstruction speed	8511 recons per second (256 ² FFT, full FoV) 40816 recons per second (256 ² FFT, 25% recFoV)
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 8 data sets

syngo MR Workplace (optional)	
Color LCD monitor and host computer as for syngo Acquisition Workplace	

In-Room syngo Acquisition Workplace (optional)	
Color LCD monitor	
Screen size (diagonal)	18"
Screen matrix	1280×1024 pixels

Euro Connector (System IEC 320)	
Load rating	100–240 Vac
	50–60 Hz
	max. 10 A

The MR system provides the level of safety according IEC 60950-1 outside patient environment. All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

Cooling system

Water cooled cabinet	Water consumption	>90 l/min ¹
	Heat dissipation to water (Maximum)	Q-engine ≤45 kW
		SQ-engine ≤60 kW
	Ready for measurement ¹	≤11.5 kW
	Typical measurement ¹	Q-engine ≤15.5 kW
		SQ-engine ≤16 kW

1. Water temperature: 6–12 °C (43–54 °F)

Line power supply

Variable	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50 Hz, 60 Hz	±1 Hz
Connection voltage ¹	Q-engine: 70 kVA SQ-engine: 110 kVA	
On-site protection	at least 160 A	

1. Long-term average in full operation

For all products, line power has to be supplied via an on-site system contactor or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710¹.

1. In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plate:

380V~ - 480V~	50/60Hz
167A - 132A	SHORT-TIME CURRENT
38A - 30A	LONG-TIME CURRENT

System off ¹	≤5.0 kW
Stand-by ¹	≤5.7 kW
Ready for measurement ¹	≤11.5 kW
Typical measurement ¹	Q-engine: ≤15.5 kW
	SQ-engine: ≤16.0 kW

1. All data incl. cold head compressor, without chilling

Dimensions

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipation [kW]
Examination Room					
Magnet 1.5 Tesla AS (incl. Helium)	207	150	210	4050±80	
Magnet in operation, incl. gradient coil, body coil, patient table, and covers	229.5	405	230	5850	3
Patient table	54 (top)	243	47–89	350	
Required min. room height clearance			235		
Min. transport dimensions	230/207 ¹	305/170 ¹	210 ¹		
Control Room					
syngo Acquisition Workplace (table+monitor)	120	80	117 (72+45)		2
Host computer	22	46	47		
syngo MR Workplace (table+monitor)	120	80	117 (72+45)		
Equipment Room					
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	160	65	189 ²	1500	5
Cooling system	65	65	189 ²	500	

1. without covers, electronics, table and table bracket

2. without attachments

- syngo MR E11 technical specifications

Magnet system

Magnet parameters	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Number of field generating coils	7
Magnet length	150 cm
Magnet bore diameter ¹	90 cm
Inner tunnel length	88 cm
Length of flared opening	54.5 cm
Inner diameter ²	60 cm
System length	160 cm
Type of installation	Fixed
Decay characteristics from full field to 20 mT	approx. 20 seconds

1. without shim coils, gradient coil, RF body coil

2. incl. shim coils, gradient coil, RF body coil

Homogeneity (based on highly accurate 24 plane plot)		
DSV	Guaranteed	Typical
50 cm	<1.5 ppm	0.8 ppm
45 cm	<1 ppm	0.4 ppm
40 cm	<0.75 ppm	0.2 ppm
30 cm	<0.5 ppm	0.1 ppm
20 cm	<0.25 ppm	0.04 ppm
10 cm	<0.05 ppm	0.01 ppm

DSV = Diameter spherical volume (x, y and z direction)

Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot meth (20 points per plane).

Standard active shim with 3 linear channels

3D Shim	
• Patient-specific automated shim • For optimization of the magnetic field homogeneity	
Time to shim	approx. 20 s

Shielding		
Active Shielding (AS)	5 th generation active shielding (AS) technology with counter coils	
Fringe field (axial×radial)	0.5 mT ¹	4.0 m×2.5 m
	0.1 mT	6.0 m×3.5 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

1. pacemaker safety limit

Magnet cooling system	
Refill interval (typical)	Not applicable ¹
Boil-off rate (typical)	0.0 l/year ¹
Helium capacity	approx. 1600 liters
Minimum helium level	40%
Cryostat	Stainless steel

1. For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

Patient handling

Tim table	
Max. patient weight (including vertical movement)	250 kg (550 lbs)
Max. scan range	154 cm optional 205 cm ¹
Vertical table movement	
Range	47–89 cm
Max. speed	10 cm/s
Horizontal table movement	
Max. range	263.5 cm
Max. speed	20 cm/s
Horizontal accuracy for repositioning from one direction	±0.8 mm

1. with Tim Whole Body Suite option

Patient comfort	
Inner diameter of magnet bore ¹	60 cm
In-bore lighting	can be set at 3 different levels
In-bore ventilation	can be set at 3 different levels
In-bore intercom	including loudspeaker, microphone and earphones

1. incl. shim coils, gradient coil, RF body coil

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2006 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise	
Patients require hearing protection with an SNR ¹ = 11 dB or more	

1. SNR = Single Number Rating

Noise for personnel in the examination room	
Noise measured	75.0 dB (A)

General: actively shielded (AS) whole-body gradient coil system

Gradient amplifier		
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses		
Q-engine	Max. output voltage ¹	1250 V
	Max. output current ¹	460 A
SQ-engine	Max. output voltage ¹	2000 V
	Max. output current ¹	625 A

1. Values for each of the 3 gradient axes

TIm [76×18] Q-engine	
Performance for each axis	
Max. amplitude	33 mT/m
Min. rise time	264 µs from 0 to 33 mT/m
Max. slew rate	125 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s

T1m [76x18] SQ-engine	
Performance for each axis	
Max. amplitude	45 mT/m ¹
Min. rise time	200 µs from 0 to 45 mT/m
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	72 mT/m
Max. eff. slew rate	346 T/m/s

1. Value is valid for longitudinal direction. Max. amplitude for horizontal and vertical direction: 40 mT/m

T1m [76x32] SQ-engine	
Performance for each axis	
Max. amplitude	45 mT/m ¹
Min. rise time	200 µs from 0 to 45 mT/m
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	72 mT/m
Max. eff. slew rate	346 T/m/s

1. Value is valid for longitudinal direction. Max. amplitude for horizontal and vertical direction: 40 mT/m

RF system

Digital signal processing system	
• Solid-state amplifier, maintenance-free • Digital quadrature demodulation • Digital filtering • Receiver with high dynamic range without adjustments	

RF transmit technology		
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.015 Hz)	
Phase control	16 bits (0.006 degrees)	
Transmitter path	Transmit amplitude	16 bit control 50 ns resolution
	Gain stability (5 min)	<0.1 dB (low signal path) <0.4 dB (power amplifier)
Transmitter amplifier	Extremely compact water-cooled	
	Transmit amplifier bandwidth	500 kHz
	Peak power	15 kW

Receive path	
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bits
Maximum number of channels ¹	76
Number of independent receiver channels that can be used simultaneously in one single scan in one FOV, each generating an independent partial image	18/32
Sampling rate (Sampling resolution)	10 MHz (100 ns)
Preamplifier noise figure	<0.5 dB
Total gain	36 dB/51 dB with automatic control
Dynamic range (referred to 1 Hz resolution bandwidth)	162 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

1. Channels (coil elements) that can be connected simultaneously

Host computer (syngo Acquisition Workplace)	
Processor	1×Intel® Xeon® Quad-Core
Clock rate	≥3.6 GHz
Main memory (RAM)	32 GB
1 st hard disk (system SW)	300 GB
2 nd hard disk (data base)	300 GB
3 rd hard disk (images)	300 GB (approx. 110000 images, matrix 256 ² and 512 ² uncompressed)
CD-R writer	approx. 4000 images 256 ² ; DICOM Standard, ISO 9660
DVD-R writer	approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
Media drives	CD/DVD drive
Color LCD monitor ¹	Screen matrix 1280×1024 pixels

1. A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

Image Processor	
Processor	≥1 Intel® Xeon® Quad-Core
Clock rate	≥3.6 GHz
Main memory (RAM)	≥16 GB
1 hard disk (system SW)	240 GB SSD
4 hard disks (raw data)	240 GB SSD
Reconstruction speed	≥8511 recons per second (256 ² FFT, full FoV) ≥40816 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 8 data sets
Media drives	CD/DVD drive

syngo MR Workplace (optional)	
Color LCD monitor and host computer as for syngo Acquisition Workplace	

In-Room syngo Acquisition Workplace (optional)	
Color LCD monitor ¹	
Screen size (diagonal)	18"
Screen matrix	1280×1024 pixels

1. A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

Euro Connector (System IEC 320)	
Load rating	100–240 Vac
	50–60 Hz
	max. 10 A

Cooling system

Water cooled cabinet	Water consumption	>90 l/min ¹
	Heat dissipation to water (Maximum)	Q-engine ≤45 kW
		SQ-engine ≤60 kW
	Ready for measurement ¹	≤11.5 kW
	Typical measurement ¹	Q-engine ≤15.5 kW
		SQ-engine ≤16 kW

1. Water temperature: 6–12 °C (43–54 °F)

Line power supply

Variable	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50 Hz, 60 Hz	±1 Hz
Connection value ¹	Q-engine: 70 kVA SQ-engine: 100 kVA	
On-site protection	at least 160 A	

1. Long-term average in full operation

For all products, line power has to be supplied via an on-site system contactor or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710¹.

1. In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plate

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT – TIME	kVA: 105
LONG – TIME	kVA: 23

Power consumption

System off ¹	≤5.0 kW
Stand-by ¹	≤5.7 kW
Ready for measurement ¹	≤11.5 kW
Typical measurement ¹	Q-engine: ≤15.5 kW
	SQ-engine: ≤16.0 kW

1. All data incl. cold head compressor, without chilling

Dimensions

The dimensions of the MR system are summarized in the following table:

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipation [kW]
Examination Room					
Magnet 1.5 Tesla AS (incl. Helium)	207	150	210	4050±80	
Magnet in operation, incl. gradient coil, body coil, patient table, and covers	229.5	405	230	5850	3
Patient table	54 (top)	243	47–89	350	
Required min. room height clearance			235		
Min. transport dimensions	230/207 ¹	305/170 ¹	210 ¹		
Control Room					
syngo Acquisition Workplace (table+monitor)	120	80	117 (72+45)		2
Host computer	22	46	47		
syngo MR Workplace (table+monitor)	120	80	117 (72+45)		
Equipment Room					
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	160	65	189 ²	1500	5
Cooling system	65	65	189 ²	500	

1. without covers, electronics, table and table bracket

2. without attachments

• SW versions differences:

Specifications	MAGNETOM Avanto (SW name and version: syngo MR B19)	MAGNETOM Avanto (SW name and version: syngo MR D13)	MAGNETOM Avanto (SW name and version: syngo MR E11)
Technical			
Software specific technologies	-	Dot (Daily Optimized Throughput) engines support	Dot (Daily Optimized Throughput) engines support
Gradient systems	Q Gradient System and 8 Channel imaging system or Q Gradient System and 18 Channel imaging system or SQ Gradient System and 18 Channel imaging system or SQ Gradient System and 32 Channel imaging system-	Q Gradient System and 18 Channel imaging system or SQ Gradient System and 18 Channel imaging system. SQ Gradient System and 32 Channel imaging system.	Q Gradient System and 18 Channel imaging system or SQ Gradient System and 18 Channel imaging system. SQ Gradient System and 32 Channel imaging system.
Operating field strength	1,5 Tesla	1,5 Tesla	1,5 Tesla
Magnet type	Superconductor	Superconductor	Superconductor
Field stability over time	<0.1 ppm/h	<0.1 ppm/h	<0.1 ppm/h
Magnet length	150 cm	150 cm	150 cm
Inner diameter	60 cm	60 cm	60 cm
System length	160 cm	160 cm	160 cm
Helium capacity	approx. 1600 liters	approx. 1600 liters	approx. 1600 liters
Minimum helium level	40%		40%
Peak power of transmitter amplifier	15 kW	15 kW	15 kW
Max. applied RF field B1+ (Body coil)	30.1 μ T	30.1 μ T	30.1 μ T
Max. specified B1+ rms (Body coil)	6.7 μ T	6.7 μ T	6.7 μ T
Max. applied RF field B1+ (Extremity coil)	47 μ T	47 μ T	47 μ T
Max. specified B1+ rms (Body coil)	17 μ T	17 μ T	17 μ T
Patient table			
Max. patient weight	250 kg	250 kg	250 kg
Max. scan range	154 cm optional 205 cm	154 cm optional 205 cm	154 cm optional 205 cm
Max. speed (vertical)	100 mm/s	100 mm/s	100 mm/s
Max. range (vertical)	470 – 890 mm	470 – 890 mm	470 – 890 mm
Max. speed (horizontal)	200 mm/s	200 mm/s	200 mm/s
Max. range (horizontal)	2635 mm	2635 mm	2635 mm
Computer			
Processor	1xIntel® Xeon® Quad-Core	Intel Xeon $\geq$ W3520 Quad-Core	1xIntel® Xeon® Quad-Core
Clock rate	$\geq$ 3.5 GHz	$\geq$ 2.66 GHz	$\geq$ 3.6 GHz
Main memory (RAM)	8 GB	$\geq$ 6 GB	32 GB

5.5 List of Coils

Coil name	Coil type	Coil design	Intended use	
			Resonant nuclei	Applications
Head Matrix	Receive	12 Channel phased array	^1H	• Head examinations • MR Angiography • Combined head/neck examinations • TMJ (temporo mandibular joints)
Neck Matrix	Receive	4 Channel phased array	^1H	• Cervical spine • Neck • Larynx/Esophagus • MR Angiography • Mediastinum • Combined head/neck examinations
Spine Matrix	Receive	24 Channel phased array	^1H	• High resolution imaging of the whole spine • Various applications in combination with additional coils
Body Matrix	Receive	6 Channel phased array	^1H	• Thorax (incl. heart) • Abdomen • Pelvis • Hip
PA Matrix (Peripheral Angio Matrix)	Receive	8 Channel phased array	^1H	• High resolution angiography of both legs incl. pelvis with high signal-to-noise ratio • Visualization of the iliac arteries and aorta • Bilateral examinations of long bones of the legs

Coil name	Coil type	Coil design	Intended use	
			Resonant nuclei	Applications
Tx/Rx CP Head Coil	Transmit/Receive	Quadrature	^1H	• Head examinations • High-resolution brain spectroscopy
32-Channel Head Array	Receive	32 Channel phased array	^1H	• High resolution head proton imaging • MR angiography of the head • Functional imaging of the brain
CP Head Array	Receive	Quadrature	^1H	• Head examinations • Head examinations with utilization of a stereotactic frame • Pediatric imaging¹
Double Loop Array	Receive	2 Channel phased array	^1H	• Simultaneous visualization of the left and right TMJs (temporo mandibular joints) with high resolution • High resolution imaging of the eye • High resolution imaging of the wrist
16-Channel Body Array anterior 16-Channel Body Array posterior	Receive	16 Channel phased array	^1H	• Thorax (incl. heart) • Abdomen • Pelvis • Hip
CP Extremity Coil	Transmit/Receive	Quadrature	^1H	• Knee • Ankle • Peripheral MR Angiography

Coil name	Coil type	Coil design	Intended use	
			Resonant nuclei	Applications
Tx/Rx 15-Channel Knee Coil	Transmit/Receive	Transmit: quadrature Receive: 15 Channel phased array	^1H	- Examinations of joints in the area of the lower extremities - High resolution knee imaging
8-Channel Knee Array	Receive	8 Channel phased array	^1H	High resolution knee or elbow imaging
Shoulder Array	Receive	4 Channel phased array	^1H	- Good visualization of small anatomical structures (e.g. labrum) - High SNR and good field homogeneity
8-Channel Wrist Coil	Receive	8 Channel phased array	^1H	High resolution wrist imaging
8-Channel Foot-Ankle Coil	Receive	8 Channel phased array	^1H	High resolution foot and ankle imaging
Breast Matrix	Receive	4 Channel phased array	^1H	- Simultaneous high resolution imaging of both breasts in all directions - Uni- or bi-lateral imaging of the breasts in sagittal direction - Visualization of axilla region
4-Channel BI Breast Coil	Receive	4 Channel phased array	^1H	- Simultaneous basic imaging of both breasts in all directions - Uni- or bi-lateral basic imaging of the breasts in sagittal direction - Uni-lateral biopsy imaging for lateral, medial and cranio-caudal access

Coil name	Coil type	Coil design	Intended use	
			Resonant nuclei	Applications
16-Channel AI Breast Coil	Receive	16 Channel phased array	^1H	• Simultaneous imaging of both breasts in all directions • Uni- or bi-lateral imaging of the breasts in sagittal direction • High resolution 2D and 3D breast imaging
2-/4-/8-Channel Sentinelle Breast Coil	Receive	Optional 2, 4, or 8 Channel phased array	^1H	• Simultaneous imaging of both breasts in all directions • Uni- or bi-lateral imaging of the breasts in sagittal direction • Uni-lateral biopsy imaging for lateral and medial access • High-resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted
2-/10-/16-Channel Sentinelle Breast Coil	Receive	Optional 2, 10, or 16 Channel phased array	^1H	• Simultaneous imaging of both breasts in all directions • Uni- or bi-lateral imaging of the breasts in sagittal direction • Uni-lateral biopsy imaging for lateral and medial access • High-resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted

Coil name	Coil type	Coil design	Intended use	
			Resonant nuclei	Applications
Prostate 2	Receive	2 Channel phased array	^1H	• Visualization of the prostate, colon, and surrounding tissue in the pelvis • Non-invasive preoperative diagnostic evaluation and treatment planning
Endorectal Coil	Receive	Linear	^1H	• Visualization of the prostate, colon, rectum, and cervix • Non-invasive preoperative diagnostic evaluation and treatment planning
4-Channel Flex Coil, large	Receive	4 Channel phased array	^1H	• Imaging of large regions such as medium to large shoulders, hip and knee
4-Channel Flex Coil, small	Receive	4 Channel phased array	^1H	• Imaging of small regions such as small to medium shoulders, wrist, elbow and ankle
Loop 11	Receive	Linear	^1H	• Examination of upper or lower extremities (e.g. shoulder, axilla)
Loop 7	Receive	Linear	^1H	• Examination of inner ear, structure of wrist and fingers, pediatrics examinations¹
Loop 4	Receive	Linear	^1H	• Examination of small structures near the surface (e.g. joints of fingers and toes, wrist, skin, temporo mandibular joints (TMJ))

Coil name	Coil type	Coil design	Intended use	
			Resonant nuclei	Applications
Flex Loop, large	Receive	Linear	^1H	- Examination of upper or lower extremities (e.g. shoulder, axilla) - Breast biopsy in combination with the optional Breast Biopsy Device
Flex Loop, small	Receive	Linear	^1H	Examination of small structures near the surface (e.g. joints of fingers and toes, wrist, skin, temporo mandibular joints (TMJ))
4-Channel Special Purpose Coil	Receive	4 Channel phased array	^1H	- Carotids - Examinations with small FOV - Small structures near the surface
CP Flex Coil, large	Receive	Quadrature	^1H	Imaging of large regions such as medium to large shoulders, hip and knee
CP Flex Coil, small	Receive	Quadrature	^1H	Imaging of small regions such as small to medium shoulders, wrist, elbow and ankle
$^{31}\text{P}/^1\text{H}$ Heart Liver Coil	Transmit/Receive	^{31}P : quadrature ^1H : linear	^{31}P ^1H	- Examinations of the heart and liver - MR imaging as well as MR spectroscopy applications

1. MR scanning has not been established as safe for imaging fetuses and infants under two years of age. The responsible physician must evaluate the benefit of the MRI examination in comparison to other imaging procedures.

5.6 Capabilities of the Medical device and ways of integration with other MDs

The following medical devices may be integrated with MAGNETOM Avanto:





	Legal Manufacturer	Way of integration
Patient table	TRUMPF Medizin Systeme GmbH + Co.KG Carl-Zeiss-Str. 7-9 07318 Saalfeld, Germany www.trumpfmedical.com	Component
Respiration cushion, respiratory belt	Pi-Products GmbH Heinrich-Hertz-Straße 8a 92224 Amberg, Germany www.pi-products.de	Accessories
Application cushion for PERU	Polyform GmbH & Co. KG Kunststofftechnik Braasstraße 15 DE-31737 Rinteln, Germany www.polyform.de	Accessories
Prostate 2 16 Ch Body Array A/P Shoulder Array Breast Array 8 Ch Knee 8 Ch Wrist 8 Ch Foot Ankle 2-/4-/8-Channel Sentinelle Breast Coil 2-/10-/16-Channel Sentinelle Breast Coil	Invivo Corporation 3545 SW 47th Ave Gainesville, Florida 32608, USA Authorized Representative (MDD93/42/EEC) GBM Authorized Representative Ltd. The White House 2 Meadow Godalming GU7 3HN, United Kingdom www.mba-gbm.com	Coil
TxRx 15 Ch Knee / TxRx CP Extremity TxRx CP Head	Quality Electrodynamics, LLC 6655 Beta Drive Suite 100, Mayfield Village, OH 44143, USA https://qedinnovations.com Authorized Representative (MDD 93/42/EEC) Medical Device Safety Service GmbH (MDSS) Schiffgraben 41 30175 Hannover, Germany	Coil
Endorectal	Bayer Medical Care, Inc. 1 Bayer Drive Indianola, PA 15051, USA Authorized Representative (MDD 93/42/EEC) Bayer Medical Care BV Horsterweg 24 6199 AC Maastricht Airport, Netherlands http://pharma.bayer.com	Coil
4 Ch BI Breast Coil	NORAS MRI products GmbH Leibnizstrasse 4 D-97204 Höchberg, Germany www.noras.de	Coil

16 Ch AI Breast Coil	RAPID Biomedical GmbH Kettelerstrasse 3-11 D-97222 Rimpfing, Germany www.rapidbiomed.com	Coil
Positioning aids (coil-specific)	Polyform GmbH & Co. KG Kunststofftechnik Braasstraße 15 DE-31737 Rinteln, Germany www.polyform.de	Accessories
Sandbags, RF homogeneity cushion, vacuum cushion, belts, table pads	Pi-Products GmbH Heinrich-Hertz Straße 8a 92224 Amberg, Germany www.pi-products.de	Accessories

5.7 Labeling information and requirements.

The following labels can be found at the electronics cabinet:

- System label

MAGNETOM Avanto (all product identifications)	
Until mid of July 2016 : (also for refurbished systems)	Starting mid of July 2016:
MAGNETOM Avanto (SW syngo MR B19)  SIEMENS MAGNETOM Avanto (1P) Model No. 07391167 (S) Serial No. 63083 Siemens AG, Wittelsbacherplatz 2, DE - 80333 Muenchen Made in Germany	MAGNETOM Avanto (SW syngo MR B19)  SIEMENS MAGNETOM Avanto (01) 04056880000011 GTIN (240) 07391167 MODEL (422) 276 ORIGIN DE (21) 125001 SERIAL (20) 119 REV E8 2015 0123 SYSTEM IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany
MAGNETOM Avanto (SW syngo MR D13)  SIEMENS MAGNETOM Avanto (1P) Model No. 10684331 (S) Serial No. 57001 Siemens AG, Wittelsbacherplatz 2, DE - 80333 Muenchen Made in Germany	MAGNETOM Avanto (SW syngo MR D13)  SIEMENS MAGNETOM Avanto (01) 04056880000026 GTIN (240) 10684331 MODEL (422) 276 ORIGIN DE (21) 157001 SERIAL (20) 12 REV 2015 0123 SYSTEM IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany
MAGNETOM Avanto (SW syngo MR E11)  SIEMENS MAGNETOM Avanto (1P) Model No. 10684332 (S) Serial No. 59001 Siemens AG, Wittelsbacherplatz 2, DE - 80333 Muenchen Made in Germany	MAGNETOM Avanto (SW syngo MR E11)  SIEMENS MAGNETOM Avanto (01) 04056880000035 GTIN (240) 10684332 MODEL (422) 276 ORIGIN DE (21) 159001 SERIAL (20) 02 REV 2015 0123 SYSTEM IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany

- Labels according to EN 50419:2006

Within the EU, products identified with this symbol are subject to guidelines 2002/96/EC for old electrical or electronics system and 2012/19/EU respectively. Please contact Siemens Service in case of questions about returning and disposing the MR system and/or its components and accessories.

- Approval identification for Canada (CSA = Canadian Standards Association)

- Approval identification for USA (UL = Underwriters Laboratories Inc.)

- Revision label

Revision:														Model - No.:																		
00	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65
66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98

The following labels can be found on the magnet cover:

- System label

MAGNETOM Avanto (all product identifications)	
Until mid of July 2016 : (also for refurbished systems)	Starting mid of July 2016:
MAGNETOM Avanto (SW syngo MR B19)	MAGNETOM Avanto (SW syngo MR B19)

MAGNETOM Avanto (all product identifications)	
Until mid of July 2016 : (also for refurbished systems)	Starting mid of July 2016:
SIEMENS MAGNETOM Avanto (1P) Model No. 07391167 (S) Serial No. 63083 Siemens AG, Wittelsbacherplatz 2, DE - 80333 Muenchen Made in Germany CE 0123 System IVK 2012	SIEMENS MAGNETOM Avanto (01) 04050000000011 GTIN (240) 07391167 MODEL (422) 276 ORIGIN DE (21) 125001 SERIAL (20) 119 REV E8 2015 0123 SYSTEM IVK Siemens Healthineers GmbH, Henkestr. 127 81062 Erlangen, Germany Made in Germany
MAGNETOM Avanto (SW syngo MR D13) SIEMENS MAGNETOM Avanto (1P) Model No. 10684331 (S) Serial No. 57001 Siemens AG, Wittelsbacherplatz 2, DE - 80333 Muenchen Made in Germany CE 0123 System IVK 2008	MAGNETOM Avanto (SW syngo MR D13) SIEMENS MAGNETOM Avanto (01) 04050000000026 GTIN (240) 10684331 MODEL (422) 276 ORIGIN DE (21) 157001 SERIAL (20) 12 REV 2015 0123 SYSTEM IVK Siemens Healthineers GmbH, Henkestr. 127 81062 Erlangen, Germany Made in Germany
MAGNETOM Avanto (SW syngo MR E11) SIEMENS MAGNETOM Avanto (1P) Model No. 10684332 (S) Serial No. 59001 Siemens AG, Wittelsbacherplatz 2, DE - 80333 Muenchen Made in Germany CE 0123 System IVK 2008	MAGNETOM Avanto (SW syngo MR E11) SIEMENS MAGNETOM Avanto (01) 04050000000035 GTIN (240) 10684332 MODEL (422) 276 ORIGIN DE (21) 159001 SERIAL (20) 02 REV 2015 0123 SYSTEM IVK Siemens Healthineers GmbH, Henkestr. 127 81062 Erlangen, Germany Made in Germany

- Approval identification for Canada (CSA = Canadian Standards Association)

- Approval identification for USA (UL = Underwriters Laboratories Inc.)

The following labels can be found at the magnet bore (front):

- Risk of injury

- Observe operator manual

The following labels can be found for laser light localizer:

- Safety certificate label (International)

- Safety certificate label (U.S.A. only)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.
Manufactured:
Siemens Aktiengesellschaft
Wittelsbacherplatz 2, D-München
Germany

Complies with 21 CFR 1040.10 and 1040.11
except for deviations pursuant to Laser
Notice No. 50, dated June 24, 2007
Manufactured: MAY 2015
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

- EAC label – for deliveries to Eurasian Economic Union

- Artwork of the MD labeling in Russian

*Siemens AG (as the Manufacturer) may be indicated on the original labeling of the Manufacturer for the systems which were produced before middle of July 2016 (also for refurbished systems);

**The labeling in Russian language is applied for the systems which will be delivered not earlier than May 2020.

For MAGNETOM Avanto with syngo MR B19 and syngo MR E11 versions:

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~3/PE 380/400/420/440/460/480 В

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6. Environment Protection Requirements.

Siemens Healthcare GmbH has a certified environmental management system, as meeting the requirements of ISO 14001:2004 for the research, development, manufacture, assembly, sale, supply, service, administration, training, serviceability, documentation and refurbishment of diagnostic imaging systems.

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8. The list of International Standards / Directives / Documents applicable to the Medical Device (ISO, CE, DoC, IEC).

EN 60601-1:2006 IEC 60601-1:2005	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
EN 60601-2-33:2011 IEC 60601-2-33:2010	Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
EN 60601-1-2:2007 IEC 60601-1-2:2007	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
EN 60601-1-6:2010 IEC 60601-1-6:2010	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
EN 60601-1-8:2007 IEC 60601-1-8:2006	Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
EN 60825-1:2014 IEC 60825-1:2014	Safety of laser products - Part 1: Equipment classification, requirements and user's guide
EN 62304:2006 IEC 62304:2006	Medical device software - Software lifecycle processes
EN 62366:2008 IEC 62366:2007	Medical devices - Application of usability engineering to medical devices
EN 1041:2008/A1:2013	Information supplied by the manufacturer of medical devices
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
EN ISO 10993- 1:2009+AC2010 ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

EN ISO 14971:2012 ISO 14971:2007	Medical devices - Application of risk management to medical devices
EN ISO 13485:2012+AC:2012	Medical devices - Quality management systems - Requirements for regulatory purposes

9. Transportation requirements.

9.1. Transportation requirements.

- The magnet can be shipped and transported worldwide with or without liquid helium. Maximum possible helium level during shipping: 95%.
- The magnet is classified as "Dangerous goods" as defined by the ADR, IMDG and IATA codes and their requirements. Liquid helium is classified under UN1963 (Helium, refrigerated liquid).
- Transportation of the system containing liquid helium must be in accordance with the appropriate national and international codes.
- Using an air-ride truck with or without a shipping pallet the following maximum peak acceleration measured at the magnet mounting points during transportation shall not exceed:
 - 5.0g vertical shock load applied to the cryostat mounting points.
 - 2.0g horizontal shock load applied to the cryostat mounting points in any direction.
 - Shock is defined as a dynamic force with duration of no more than 0.5 seconds.
- Road transport is allowed up to 15.000 miles.
- Road transport may be with or without an approved shipping pallet. Without NO shipping pallet, the magnet must be supported using transport beams.
- Road vehicles must be fitted with air ride suspension, whether or not a shipping pallet is used.
- Always transport a magnet using vehicles fitted with air ride suspension.
- Road vehicles must be suitable for the intended usage, be properly maintained at all times, and comply with all local laws and regulations.
- Relevant requirements of the ADR code shall be satisfied.
- The vehicle and operator must be approved by Siemens.
- This specification does not constitute acceptance for use as part of a mobile system. Mobile magnets require a specific trailer design for continuous road transportation.
- Air transport is allowed up to 30.000 miles.
- If air transport is used for any section of the route, an approved shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IATA code must be satisfied.
- Sea transport is allowed up to 30.000 miles.

- If sea shipment is used for any section of the route, an approved (magnet system related) shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IMDG code must be satisfied.
- Rail Transportation is not allowed!
- Transportation through the channel tunnel (UK to France) is permissible if the system is mounted on an approved air ride truck. The air ride suspension may be inactive.
- Lifting the magnet must always be conducted by personnel who are trained in handling components and medical equipment of this size and weight.
- Lifting during transportation including the final movement onto the installed site must be performed using a suitable sized crane capable of lifting the full magnet weight.
- Suitable straps or chains capable of lifting the full magnet weight must be used. The straps or chains must be secured to the magnet lifting points using the four shackles supplied and fitted to the magnet.
- All transport activities including delivery of the magnet and the system components to the MR suite may be performed by a licensed carrier only.

9.2. Transportation specifications – Magnet

Transport specifications for the Magnet	
Acceleration amplitude	2 m/s ² (2 - 100 Hz)
Peak acceleration	5 g vertical, 2 g horizontal (20 ms time period)
Storage temperature	-25 °C to +60 °C
Relative humidity	10% to 95% non-condensing
Atmospheric pressure	20 kPa to 106 kPa

9.3. Transportation specifications – System Components

Transport specifications for the system components	
Acceleration amplitude	2 m/s ² (2 - 100 Hz)
Peak acceleration	5 g vertical, 2 g horizontal (20 ms time period)
Ambient temperature	-25 °C to +60 °C
Relative humidity	10% to 95%
Atmospheric pressure	70 kPa to 105 kPa

10. Storage and Service Life.

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order, and is therefore not susceptible to degradation that would lead to functional failure.

Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR scanner. The MR scanners do not spend time on the “shelf” that would pose a risk to the performance of the device or its components.

The expected service life of 10 years of operation under the environmental conditions indicated in the instructions for use as well as with all preventive and corrective maintenance performed by qualified technical personnel authorized by Siemens Healthcare GmbH.

Refurbished systems are considered not to be different compared to standard (non-refurbished) systems regarding the expected service life, since production- and test-procedures including service and maintenance are identical to systems which have not been refurbished.

If the systems are used beyond this time, additional maintenance and service may be required in order to maintain system performance and reliability.

10.1 Storage specifications – Magnet

Storage specifications for the Magnet	
Acceleration amplitude	2 m/s ² (2 - 100 Hz)
Peak acceleration	5 g vertical, 2 g horizontal (20 ms time period)
Storage temperature	-20 °C to +60 °C
Relative humidity	10% to 95% non-condensing
Atmospheric pressure	20 kPa to 106 kPa

10.2 Storage specifications – System Components

Storage specifications for the system components	
Acceleration amplitude	2 m/s ² (2 - 100 Hz)
Peak acceleration	5 g vertical, 2 g horizontal (20 ms time period)
Ambient temperature	-25 °C to +60 °C
Relative humidity	10% to 95%
Atmospheric pressure	70 kPa to 105 kPa

11. Procedure of recycling and disposal (Medical Device recycling or disposal data).

On disposing of the complete system or parts thereof, currently valid environmental legislations must be observed. Examples of environmentally relevant components are:

- Accumulators and batteries;
- Transformers;
- Capacitors.

All medical waste shall be collected, used, decontaminated, placed, stored, transported, listed and disposed of in accordance with the laws on public health and epidemiological welfare.

System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

For details the local customer service representative shall be contacted.

12. Usage conditions.

12.1. Installation requirements.

To ensure correct installation, complete that all preliminary steps specified. For detailed information on preparing the room, transporting the system, and other installation requirements, consult the project manager. The on-site project manager is responsible for all room and system-related preparation.

The following prerequisites must be fulfilled:

- Ensure that the magnet room is properly equipped.
- Ensure that the openings of the RF shield are sufficiently large (according to the height and width specified in the Planning Guide).
- Ensure that there is sufficient space to install the magnet.
- Take account of narrow passages due to suspended ceilings, carriers, and columns.
- Determine the final location of the magnet in the magnet room.
- Ensure that the primary cold water supply is ready for operation.
- Block off the unloading zone to ensure that passersby are not harmed.
- Do not block escape routes or essential access routes to the hospital.
- The rooms must be properly equipped. The RF shield, including the opening for inserting the magnet, must be installed. The magnet room must also be equipped with sufficient lighting. The floor must be clear so that the magnet and system components can be transported.
- The four corners of the magnet must be drawn on the floor of the magnet room to determine the magnet's location. Refer to the project plan.
- The floor from the unloading zone to the final location must be covered with metal plates (at least 5 mm thick).
- The tools required for installation must be available on-site (steel rollers, hydraulic lifts, etc.).

12.2. Information about maintenance (general and on the territory of Russia).

Routine checks and periodic maintenance are necessary to ensure safe and satisfactory operation of the system. Proof of these activities is required by authorities in certain countries.

Routine checks include:

Daily, weekly, and monthly inspections as well as legally required checks as described in the system "Function and Safety Checks" chapter of the operating instructions. Normally, the system operator entrusts the clinical operating personnel with the task of performing these routine checks.

Periodic maintenance includes:

- Safety check (including safety-related tests)
- Preventive maintenance
- Quality and function tests
- Replacement of safety-relevant parts subject to wear.

This work has to be performed by qualified and authorized service engineers only. Qualified in this context means that the engineers have been trained accordingly or have acquired practical experience through routine service activities. Authorized means that the engineers have been authorized by the operator of the system to perform maintenance work.

Upon first start-up of the system, Siemens ask the system operator to designate a staff member responsible for ensuring that routine checks, preventive inspection, and maintenance work are performed. This staff member is responsible for archiving all certificates in the "System Owner Manual" binder.

In addition to our repair service, Siemens also offers a complete range of services for the preventive inspection and maintenance of the system. These services can be called on as required or agreed upon in a flexibly drafted maintenance contract.

Service Contacts are:

Siemens Healthcare GmbH – Customer Service (CS)
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: siemens-healthineers.com

Siemens Healthcare LLC
Address: Dubininskaya Str. 96, 115093 Moscow, Russia
Tel.: +7 495 737 12 52
Web: siemens-healthineers.com/ru/

Customer's support center of Service Department Siemens Healthcare LLC
8 800 550 50 22
Service.ru@siemens-healthineers.com

12.3. Information about calibration of the MD / requirements to calibration (if applicable).

There is no calibration required for the MR systems. Nevertheless there are different checks and tests which have to be done as described in the system owner manual.

13. Warranty

The system owner is responsible for the software licenses including software warranty conditions.

The specific warranty conditions are usually defined in the contractual agreement between Siemens Healthcare GmbH and the system owner.

14. Contact for reclamation

In case of any questions or reclamations, the following shall be contacted:

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Herewith, I certify the foregoing signature,
acknowledged before me, of

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
Mrs. Velasco Alvarado is personally
known.

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