SEELER & DR.FELLMETH
NOTARE

NOTARE SEELER & DR. FELLMETH · Q7, 23 · 68161 MANNHEIM

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D-Dimer Gen. 2 Calibrator Set

cobas®

- 2 Amiral J, Plassart V, Minard F. Measurement and clinical relevance of D-dimer by ELISA. In: fibrinogen and its derivatives. Müller-Berghaus G, Scheefers-Berchel U, Selmayr E and Henschen A, eds. 285-290 experta Medica, Amsterdam 1986.
- 3 Occupational Safety and Health Standards: Bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
- 4 Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.

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1. Adema E, Gebert U. Pooled patient samples as reference material for D-Dimer. *Thromb Res* 1995;80(1):85-88.
2. Amiral J, Plassart V, Minard F. Measurement and clinical relevance of D-dimer by ELISA. In: *fibrinogen and its derivatives*. Müller-Berghaus G, Scheefers-Berchel U, Selmayr E and Henschen A, eds. 285-290 *experta Medica*, Amsterdam 1986.
3. Occupational Safety and Health Standards: Bloodborne pathogens. (29 CFR Part 1910.1030). *Fed. Register*.
4. Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.

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Design and Manufacturing

Active Ingredients, Biological Safety

D-DI2 (**REF** 07429410190)

D-DI2 Cal Set (**REF** 07683456190)

D-DI2 Con (**REF** 07571933190)

1) *Active Ingredients*

List of active ingredients (see certificates attached):

- D-DI2**
 - Latex particles coated with monoclonal anti-human D-Dimer antibodies (mouse)
- D-DI2 Cal Set**
 - Human Serum with human D-Dimer fragments
- D-DI2 Con**
 - Human Serum with human D-Dimer fragments

2) *Biological Safety*

General statement regarding biological risk of IVDs

An In vitro diagnostic (IVD) medical device is intended by the manufacturer to be used in vitro for the examination of specimens, including blood and tissue donations, derived from the human body.

The major risk in the use of IVDs is associated with patient samples which may contain infectious agents.

CONFIDENTIAL AND PROPRIETARY

Doc.No.: 07429410190-07683456190-07571933190_DM-Actingr_DM-Biolog_12-Jun-2017

In order to minimize this risk, proper handling is recommended in all instructions for use under the section *Precautions and warnings*.

For manufacturing of D-DI2, D-DI2 CI Set and D-DI2 Con no material of recombinant origin is used.

Elimination or prevention of risk for IVD users related to biological material

Materials of Human origin

All materials of human origin used for the manufacturing of IVDs (mostly manufacture of calibrators and controls) have to be tested negative/non-reactive as described in the respective IFU. Testing is done inhouse or by the supplier (see certificates attached)

Nevertheless, all blood-derived materials should be handled as potentially infectious.

Materials of animal origin

Certificates by the suppliers are provided (see certificates attached).

Biological Material used

Please refer to the composition.

CONFIDENTIAL AND PROPRIETARY

Doc.No.: 07429410190-07683456190-07571933190_DM-Actlgr_DM-Biolog_12-Jun-2017

Certificate of Analysis

We hereby certify that the following active ingredient, which is used for D-Dimer R2 bulk reagent, met our specification.

Active ingredient: DD10075J21010 Anti-D-Dimer antibody

Lot no. 20140301

July 08, 2014

Parameter	Results
1. Appearance:	yes
Method: sensory evaluation	
Judgment: White latex suspension	
2. Accuracy:	
Method: Assaying the QC samples	
Judgment:	
1) within +/-30% of values assayed by the previous lot	QC-1 +0.01 μ gFEU/mL
	QC-2 101 %
	QC-3 102 %
	QC-4 103 %
	QC-5 102 %
2) within +/-10% of the target values > 1.7 μ gFEU/mL	QC-1 +0.13 μ gFEU/mL
+/-0.17 μ gFEU/mL of the target values $\leq$ 1.7 μ gFEU/mL	QC-2 105 %
	QC-3 102 %
	QC-4 100 %
	QC-5 103 %

This lot met our specification.

 **LSI Medience Corporation**

THE KAITOKI Bldg., 13-4, Uchikanda 1-chome, Chiyoda-ku
Tokyo 101-8517, Japan


THE KAITOKI COMPANY

ASSOCIATED CHEMICAL INDUSTRIES GROUP

Production date: July 07, 2014

Expiry date: August 08, 2021

LSI Medience Corporation

Yuichi Shintani, General Manager

Quality Control Department

Certificate of Analysis

We hereby certify that the following raw material, which is used for D-Dimer Gen.2 bulk reagent, met our specification.

Raw material of animal origin : Gamma Globulin from Mouse Serum

Lot no. 20150801

September 16, 2015

Item	Method	Requirement	Results
Appearance	Visual Contact	Solid (powder)	Pass
Purity	HPLC	80% or more	89% Pass

This lot met our specification.

Expiry date: August 03, 2025

LSI Medience Corporation

Shigekazu Sogabe, General Manager
Production and Logistics Center

Certificate of Analysis

Active ingredient : DD10075P01210 human D-Dimer fragment

Lot no. 20140801

November 27, 2014

Parameter	Results
1. Appearance:	yes
Method: sensory evaluation	
Judgment: transparent and colorless solution	
2. Accuracy:	no low molecular D-Dimer
Method: SDS Page	
Judgment: no low molecular D-Dimer	
3. Safety:	negative
Method: Checking the safety certificate issued by manufacturer	
Judgment: HCV,HB and HIV	

This lot meets our specification.

Production date: November 27, 2014

Expiry date: September 30, 2034

LSI Medience Corporation

Shigekazu Sogabe, General Manager
Production and Logistics Center

Certificate of Analysis

We hereby certify that the following raw material, which is used for D-Dimer Gen.2 bulk reagent, met our specification.

Raw material of animal origin: Albumin from Bovine Serum

Lot no. 20150101

January 20, 2015

Item	Method	Requirement	Results
Appearance	Visual Contact	Solid (powder)	Pass
Purity	Electrophoresis	98% or more	100% Pass

This lot met our specification.

Expiry date: June 20, 2020

LSI Medience Corporation

Shigekazu Sogabe, General Manager
Production and Logistics Center

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07683456190

Scope of the Analysis

The following areas of product risks were considered:

1 Health

1.1 Indirect hazards to patient health caused by assay results

Confidential
- For registration purposes only -

- 1.1.1 Errors in design
- 1.1.2 Errors in development (performance data)
- 1.1.3 Faulty reagent
 - 1.1.3.1 Errors during use
 - 1.1.3.2 Errors in production
- 1.1.4 Faulty calibrator
 - 1.1.4.1 Errors during use
 - 1.1.4.2 Errors in production
- 1.1.5 Faulty control
 - 1.1.5.1 Errors during use
 - 1.1.5.2 Errors in production
- 1.1.6 Sample
- 1.1.7 Errors in the results evaluation or interpretation
- 1.2 Injury to the user's health by reagents, calibrators and controls
- 1.3 Injury to manufacturer's health by reagents, calibrators and controls
- 2 Environment
 - Reagents, calibrators and controls
- 3 Property
 - Reagents, calibrators and controls

Environmental influences may apply by temperature or light during shipment, storage and application. Stability is a potential risk addressed in the risk area "Incorrect assay results".

Microbial susceptibility is investigated by routinely established procedures during the design. In addition, where indicated, microbial status is observed as part of the production monitoring.

The In-vitro diagnostic assay device is intended for application on the respective instruments as described in the instruction for use. Potential failure modes of instruments of Centralised and Point of Care Solutions and the potential consequences on test results are investigated in the corresponding analysis for the respective instrument. Potential risks to property will be addressed by this analysis too.

Specifically for D-Dimer:

Hazard	Hazardous situation	Risk Mitigation
<i>Incorrect assay result</i>	<i>Li-Heparin plasma. (suitable sample material on ClinChem analyzers) may cause wrong results on Coag analyzers by sample carry-over</i>	<i>Limit usage of sample material for D-Dimer determination on Coag analyzers to citrated plasma and include information in IFU/method sheet.</i>

Confidential
- For registration purposes only -

Hazard	Hazardous situation	Risk Mitigation
Immunoglobulin interference	Interference is more pronounced than on clin. Chem. Instruments	Pipetting scheme is identical to Clin.-Chem. instruments, therefore effectiveness of interference reduction is likely the same. Will be tested during application work if possible. Statement in IFU like for Clin.-Chem. Instruments
		Residual risk is accepted as it is state of the art for immunological assays

There have been 0 residual risks in the red area after applying risk control, that have to be investigated by a specific detailed risk - benefit analysis.

Summary Statement and Acceptance

During the analysis we investigated the potential risks along the scope described above.

By application of Internally established methods of risk control all these risks could be reduced to an acceptable value. Those residual risks that may be of relevance to the user, for example limitations in performance, special guidance for handling of the product, dangerous substances, etc. are described in the instructions for use or by means of warnings.

As a result of the analysis, it is concluded, that the product conforms to safety principles, taking into account the generally acknowledged state of the art,

- when used under the conditions and for the purpose intended
- when stored and handled taking account of the instructions and information provided by the manufacturer

and that the overall benefit of the product outweighs the residual risks.

The risk management during the product lifecycle underlies the regulations of the QM-system of Centralised and Point of Care Solutions.

Author (Project-/Development Lead or Product Responsible Person):

Dr. Sylvia Rosenblatt

Date (dd-mmm-yyyy):

06-Mar-2017

Confidential

- For registration purposes only -

Signature:

A handwritten signature in black ink, appearing to be "S. H. ...", written over a horizontal line.

If this document is released within a Document Management System no entry is required here.

Confidential
- For registration purposes only -

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Mannheim, den 12. FEB. 2021

Claudia Seeler
Notarin in Mannheim

For Presentation to the Authorities of Russia

Explanation Letter

We,
Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim
Germany,

are the legal manufacturer for the medical devices:

- «Reagent in cassette for the quantitative determination of fibrin degradation products (D-dimer and X-oligomer) by latex-enhanced immunoturbidimetry in citrated plasma on cobas t automated hemostasis analyzers (D-DI 2 / Tina-quant D-Dimer Gen.2 cobas t systems)" manufactured by Roche Diagnostics GmbH, Ref. No. 81981 dated October 19, 2020»
- «Calibrator set for the quantitative determination of fibrin degradation products (D-dimer and X-oligomer) by latex-enhanced immunoturbidimetry in citrated plasma on cobas t automated hemostasis analyzers (D-DI2 Cal Set / D-Dimer Gen. 2 Calibrator Set cobas t systems)" manufactured by Roche Diagnostics GmbH, Ref. No. 81989 dated October 19, 2020»
- «Set of control materials for the quantitative determination of fibrin degradation products (D-dimer and X-oligomer) by latex-enhanced immunoturbidimetry in citrated plasma on cobas t automated hemostasis analyzers (D-DI2 Con / D-Dimer Gen.2 Control I/II cobas t systems)" manufactured by Roche Diagnostics GmbH, Ref. No. 81985 dated October 19, 2020»

Express to you our respect and due to requests № 10-73601/20 from 11.12.2020, №10-74449/20 from 15.12.2020 and №10-74450/20 from 15/12/2020 we would like to inform you about the following:

The abovementioned products are dedicated for the quantitative immunological determination of fibrin degradation products (D-dimer and X-oligomers) in human citrated plasma on the indicated cobas t analyzer.


Fibrin clot formation:

Thrombin converts fibrinogen to soluble fibrin by cleaving the fibrinopeptides A and B. The fibrin monomers polymerize spontaneously. Active factor XIII links two D-domains and generates a solid fibrin clot. Thereby a new plasmin-resistant antigenic determinant ("D-dimer") is produced, which can be detected immunologically.

Fibrin clot degradation:

During the plasmin induced degradation process of a fibrin clot many different sized fragments containing D-dimer determinants are formed. These high molecular weight fibrin degradation products are referred to as X-oligomers. The above-mentioned assay has a strong affinity to the D-dimer determinant and therefore also detects X-oligomers.

From [1]: A variety of different fibrin degradation products containing D-Dimer determinants (D-D) are formed during digestion.

Quantitative determination of D-dimer fragments:

Latex particles are coated with monoclonal antibodies (F(ab')₂ fragments) to the D-dimer epitope. The antigen/antibody complexes produced by the addition of samples containing D-dimer determinants (which includes X-oligomers) leads to an increase in the turbidity of the test. The change of absorbance is determined over time correlates with the concentration of D-dimer epitopes in the sample.

As mentioned in [2]: "The degradation products from fibrin clot lysis are a heterogeneous group of substances, [...]. Usually the term "D-dimer" is a collective term, in reference to this entire group of fragments."

For descriptive reasons we decided to refer to the term "X-oligomers" only once in the Intended Use. For the remaining chapters of the Package Inserts the commonly used scientific collective abbreviation "D-dimer" was chosen.

Due to that in "Specification of intended use section" in Annex to Instructions for use, only D-Dimer determination could be mentioned, without information on determination of X-oligomers:

In vitro assay for the quantitative immunological determination of fibrin degradation products (D-dimer) in human citrated plasma on the indicated cobas t analyzers.

In conjunction with a non-high clinical probability assessment, a normal D-dimer result (< 0.5 µg FEU^a/mL) is intended as an aid in the exclusion of deep vein thrombosis (DVT) and pulmonary embolism (PE).

a) Fibrinogen Equivalent Unit

Finally, we ask you to use the following names of medical devices:

- «Reagent in cassette for the quantitative determination of D-dimer by latex-enhanced immunoturbidimetry in citrated plasma on cobas t automated hemostasis

analyzers (D-DI 2 / Tina-quant D-Dimer Gen.2 cobas t systems)" manufactured by Roche Diagnostics GmbH »

- «Calibrator set for the quantitative determination of D-dimer by latex-enhanced immunoturbidimetry in citrated plasma on cobas t automated hemostasis analyzers (D-DI2 Cal Set / D-Dimer Gen. 2 Calibrator Set cobas t systems)" manufactured by Roche Diagnostics GmbH»
- «Set of control materials for the quantitative determination of D-dimer by latex-enhanced immunoturbidimetry in citrated plasma on cobas t automated hemostasis analyzers (D-DI2 Con / D-Dimer Gen.2 Control I/II cobas t systems)" manufactured by Roche Diagnostics GmbH »

Roche Diagnostics GmbH is responsible for quality of the product and believes that the mentioned information minimizes the risk of user misinformation and does not entail any changes in the properties and characteristics that affect to the quality, efficacy and safety of the above mentioned medical products, produced earlier and currently. The functionality and operation of the above mentioned medical devices, produced earlier and currently, do not change.

Therefore, we ask you to continue the registration process of the medical device.

Roche Diagnostics GmbH is responsible for design and development, manufacturing, post-market surveillance and quality control of all products manufactured.

Roche Diagnostics GmbH guarantees quality, efficacy and safety of a medical product.

References:

- [1] Favresse J, Lippi G, Roy PM, Chatelain B, Jacqmin H, ten Cate H & Mullier F. D-dimer: Preanalytical, analytical, postanalytical variables, and clinical applications. Critical Reviews in Clinical Laboratory Sciences; Volume 55; Issue 8: 548-577 (2018)
- [2] Stang LJ. D-dimer and fibrinogen/fibrin degradation products. Haemostasis: methods and protocols (ed. Monagle P). Methods in molecular biology 2013;992:415-42

Mannheim, February 11, 2021

Sincerely,


Roche Diagnostics GmbH

i.V. / on behalf of the company

i.V. / on behalf of the company

Andrea Weber
Manager Regulatory Affairs

Dr. Manfred Böhm
Manager Regulatory Affairs

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Explanation Letter

We,
Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim
Germany,

Explanation of shelf-life determination procedure for Roche Lab Coagulation products

To whom it may concern.


The purpose of the real time stability study is to examine the possibility of long-term storage of medical devices for IVD diagnostics at declared temperature range.

The testing of the real-time stability is done in full compliance with EN ISO 23640. The function of the reagent's lot is tested at least at 4 measurement points, including the beginning and the end. The shelf life will be claimed on the basis of the results for real-time stability testing for all three lots minus a safety margin.


Shelf life is determined by testing three lots using the following model:

Shelf life determination – Model 3 “Real-time”

• Real-time with transport + “in-use” (R&D or QC)

• Real-time (R&D or QC)

• Declaration: (R&D or QC)

Stability confirmed on 3 lots (minus the safety margin)

1

Picture 1: Shelf life determination – Model 3 „Real Time”

Data from three lots in real time are available. One of the lots accounts for the transport stress and "in-use" stress in the form of a "worst case" stress condition. The shelf life to be assigned is based on the lot having the shortest demonstrated stability minus the safety margin.

The testing procedure for stability testing is the functional test. Acceptance criteria are determined according internal quality control procedure, production specifics and test constructions. The results of this real time stability program display the stability under the tested conditions. For product related information, please refer to the corresponding Stability Certificate.

For Lab Coagulation portfolio, it was decided to use shelf life testing according model 3 (picture 1). The accelerated shelf life determination mode is used to predict the stability and in any case, it must be confirmed using realtime data based on model 3.

Due to that, we ask you to continue the registration process of the medical devices with the abovementioned information.

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
Roche Diagnostics GmbH guarantees quality, efficacy and safety of medical products.

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Explanation Letter

We,
Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim
Germany,

Declaration of Minimum Shelf-life (MSL) of Roche Lab Coagulation products

To whom it may concern.

Roche Diagnostics GmbH, Sandhofer Strasse 116, 68305 Mannheim, Germany hereby inform that the shelflife of all Roche products, of which we are the legal manufacturer, is an indication of the date until when the IVD medical device may be used safely, expressed as the year, month, and day.


The standard manufacturing process employs the following activities before the products are ready for sale and distribution to the countries:

- 1) Assignment of calibration-curves and calibration-values
- 2) In-process control according to pertinent testing procedure
- 3) Release-testing of intermediated according to pertinent testing procedure
- 4) Functional testing of assay-performance according to pertinent testing procedure
- 5) Packaging and labelling of product
- 6) Batch-Record-Review and final release of product to the market by Quality Assurance

The above mentioned activities may take several months from manufacturing date to distribution. Hence, a product manufactured on a particular date is usually available for dispatch/batch release after a time lag of several months. The respective Certificate of Analysis (CoA) of the pertinent product reveals this time lag, as the span between the Date of Manufacture and the Date of Batch-Release i.e. "Issue Date", also imprinted in the CoA.

By virtue of the final release, a product is only available to the market after it has been approved and its Time-in-Market is a fraction of the registered (total) shelf life of the product.

The facts mentioned above are shown at the following scheme.

* Registered Shelf Life can be equal to, or longer than the lot-specific kit shelf life

- DoM (Date of Manufacturing) is specific to a kit lot
- Kit lot is the 6 or 8-digit identifier that is specific to a single kit-lot.
- Components are the individual reagents that make up a kit.
- Registered Shelf life is the maximum allowable shelf life of a kit (stated in the Stability Certificate)
- Kit Date of Expiry (DoE) is the calculated kit expiry date based on the shelf life of the component with the shortest expiry date.
- Date of Landing (DoL) is the date that a kit is received at the port of entry
- Residual Shelf Life is the period of time between Date of Landing and DoE
- Calculated shelf is the lot-specific shelf life that is calculated by DoE — DoM
- SSLED (Start of Shelf Life Expiry Date) is used as the starting point for calculating the DoE based on the registered shelf life.
- The period of time between the DoM and the DoE of a kit may be shorter than the Registered Shelf Life
- Residual shelf life is the amount of remaining shelf life for a kit lot (lot-specific DoE minus lot-specific DoM).
- Calculating the residual shelf life as a percentage of the registered shelf life:

$$\left(\frac{\text{Lot-specific DoE} - \text{DoL}^*}{\text{(in months)}} \div \frac{\text{Registered shelf life}}{\text{(in months)}} \right) \times 100 \Rightarrow \text{Residual shelf life as \% of Registered Shelf Life}$$

In the registration dossier Registered Shelf life (maximum allowable shelf life of a kit) is listed.

Considering above mentioned facts and with the intent to make the products available with suitable date of expiration on the labeling the lag between Date of Manufacture and Date of Expiry can be either equals to Registered Shelf life or less than it.

We hereby confirm that the Time-in-Market after accommodation of the time lag, due to activities in manufacturing as referred above, will have no detrimental impact on the quality or performance of the products.

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Roche Diagnostics GmbH guarantees quality, efficacy and safety of all medical products.


Mannheim, 10 December 2020

Sincerely,

Roche Diagnostics GmbH

i.V./on behalf of the company

i.V./on behalf of the company

Andrea Weber


Manager Global Regulatory Affairs

Christian Bovenkerk

Project Manager Regulatory Affairs

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A blue ink signature of L.V. Dejneko, written in a cursive style.

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27. Mai 2020

Signature

Andrea Weber
Manager Global Regulatory
Affairs

Critical Reviews in Clinical Laboratory Sciences

ISSN: 1040-8363 (Print) 1549-781X (Online) Journal homepage: <https://www.tandfonline.com/loi/ilab20>

D-dimer: Preanalytical, analytical, postanalytical variables, and clinical applications

Julien Favresse, Giuseppe Lippi, Pierre-Marie Roy, Bernard Chatelain,
Hugues Jacqmin, Hugo ten Cate & François Mullier

To cite this article: Julien Favresse, Giuseppe Lippi, Pierre-Marie Roy, Bernard Chatelain, Hugues Jacqmin, Hugo ten Cate & François Mullier (2018) D-dimer: Preanalytical, analytical, postanalytical variables, and clinical applications, Critical Reviews in Clinical Laboratory Sciences, 55:8, 548-577, DOI: [10.1080/10408363.2018.1529734](https://doi.org/10.1080/10408363.2018.1529734)

To link to this article: <https://doi.org/10.1080/10408363.2018.1529734>

© 2019 The Author(s). Published by Informa UK Limited, trading as Taylor & Francis Group.

Published online: 10 Jan 2019.

Submit your article to this journal [↗](#)

Article views: 3455

View related articles [↗](#)

View Crossmark data [↗](#)

Citing articles: 1 View citing articles [↗](#)

Full Terms & Conditions of access and use can be found at
<https://www.tandfonline.com/action/journalInformation?journalCode=ilab20>

D-dimer: Preanalytical, analytical, postanalytical variables, and clinical applications

Julien Favresse^a, Giuseppe Lippi^b , Pierre-Marie Roy^c, Bernard Chatelain^a, Hugues Jacqmin^a, Hugo ten Cate^d and François Mullier^a 

^aCHU UCL Namur, Namur Thrombosis and Hemostasis Center, Hematology Laboratory, Université Catholique de Louvain, Yvoir, Belgium; ^bSection of Clinical Biochemistry, University Hospital of Verona, Verona, Italy; ^cDépartement de Médecine d'Urgence, CHU d'Angers, Institut MITOVASC, Université d'Angers, Angers, France; ^dDepartment of Internal Medicine, Cardiovascular Institute, Maastricht University Medical Center, Maastricht, the Netherlands

ABSTRACT

D-dimer is a soluble fibrin degradation product deriving from the plasmin-mediated degradation of cross-linked fibrin. D-dimer can hence be considered a biomarker of activation of coagulation and fibrinolysis, and it is routinely used for ruling out venous thromboembolism (VTE). D-dimer is increasingly used to assess the risk of VTE recurrence and to help define the optimal duration of anticoagulation treatment in patients with VTE, for diagnosing disseminated intravascular coagulation, and for screening medical patients at increased risk of VTE. This review is aimed at (1) revising the definition of D-dimer; (2) discussing preanalytical variables affecting the measurement of D-dimer; (3) reviewing and comparing assay performance and some postanalytical variables (e.g. different units and age-adjusted cutoffs); and (4) discussing the use of D-dimer measurement across different clinical settings.

Abbreviations: AAD: acute aortic dissection; ACCF: American College of Cardiology Foundation; AD: aortic dissection; AHA: American Heart Association; AMI: acute mesenteric ischemia; APTT: activated partial thromboplastin time; ASA: American Stroke Association; BMI: body mass index; CAP: College of American Pathologists; CDR: clinical decision rule; CLIA: chemiluminescent enzyme immunoassay; CLSI: Clinical Laboratory Standards Institute; CPR: clinical prediction rules; CT: computed tomography; CTPA: CT pulmonary angiogram; CV: coefficient of variation; CVT: cerebral venous thrombosis; DD/E: D-dimer/fragment E complex; DDU: D-dimer units; DIC: disseminated intravascular coagulation; DiPEP: Diagnosis of PE in Pregnancy; DVT: deep venous thrombosis; eGFR: estimated glomerular filtration rate; ELFA: enzyme-linked fluorescence immunoassay; ELISA: enzyme-linked immunosorbent assays; ESC: European Society of Cardiology; ESTES: European Society for trauma and Emergency Surgery; F/T: freezing/thawing; F1 + 2: prothrombin fragments 1 + 2; FACT: fibrin assay comparison trial; FDA: Food and Drug Administration; FDP: fibrin degradation products; FEU: fibrinogen equivalent units; FV: factor V; FVIII: factor VIII; G: Gauge; GP: general practitioner; HMWF: high molecular weight fibrinogen; HR: hazard ratio; HV: healthy volunteers; ISO: International Organization for Standardization; ISTH: International Society on Thrombosis and Haemostasis; IV: intravenous; LMWF: low molecular weight fibrinogen; NNT: number needed to test; NPV: negative predictive value; PE: pulmonary embolism; PET: polyethylene terephthalate; POC: point of care; PPV: positive predictive value; PT: prothrombin time; PTS: pneumatic tube system; RS: risk score; RT: room temperature; TAT: turnaround time; t-PA: tissue plasminogen activator; TT: thrombin time; URL: upper reference limit; VKA: vitamin K antagonist; VTE: venous thromboembolism; WHO: World Health Organization

ARTICLE HISTORY

Received 23 June 2018
Revised 31 August 2018
Accepted 25 September 2018
Published online 10 January 2019

KEYWORDS

D-dimer; preanalytical;
analytical; postanalytical;
venous thromboembolism

Introduction: what is D-dimer?

D-Dimer is a generic term referring to multiple peptide fragments deriving from plasmin-mediated degradation of cross-linked fibrin. Their presence reflects concomitant activation of both coagulation and fibrinolysis

[1–4]. The first step of D-dimer formation involves the thrombin-catalyzed conversion of fibrinogen into fibrin monomers (Figure 1). Human fibrinogen is a soluble plasma glycoprotein composed of three different pairs of polypeptide chains (α -, β -, and γ -) connecting

CONTACT François Mullier  francois.mullier@uclouvain.be  CHU UCL Namur, Hematology Laboratory, Université Catholique de Louvain, Avenue G. Thérèse, 1, 5530 Yvoir, Belgium

© 2019 The Author(s). Published by Informa UK Limited, trading as Taylor & Francis Group.

This is an Open Access article distributed under the terms of the Creative Commons Attribution-NonCommercial-NoDerivatives License (<http://creativecommons.org/licenses/by-nc-nd/4.0/>), which permits non-commercial re-use, distribution, and reproduction in any medium, provided the original work is properly cited, and is not altered, transformed, or built upon in any way.


Figure 1. Mechanism of D-dimer production [1,3,4,164].

two outer D-domains to the central E-domain [5,6]. Thrombin enzymatically cleaves two cryptic polymerization sites located on the E-domain, thus leading to the generation of both highly self-adhesive fibrin monomers and fibrinopeptides A and B [1]. Fibrin monomers then bind to one another to form a soluble network [1]. Simultaneously, the complex between soluble fibrin polymers, thrombin, and plasma factor XIII promotes the formation of factor XIIIa, which catalyzes covalent cross-linking of fibrin polymer via intermolecular bonds formed between lysine and glutamine residues, thus enabling the generation of stable and insoluble clots [1,4]. Later, the fibrinolytic pathway leads to degradation of stabilized clots through plasmin activation [4]. Plasmin is activated from plasminogen by tissue plasminogen activator (t-PA) at the fibrin surface and cleaves fibrin at specific sites [7]. The products of this reaction are the conventionally-defined fibrin degradation products (FDP), and display a vast array of molecular weights. In the first steps of fibrinolysis, FDP are large. Continual breakdown generates the fragment D-dimer/fragment E complex (DD/E), which has two covalently-bound D-domains [4]. Because of variable degrees of plasmin-mediated proteolysis, plasma samples then contain a mixture of fibrin fragment complexes containing one or several D-dimer motifs whose molecular weights range from 228 kDa (DD/E) to several thousand kDa (X-oligomers) [8,9]. These FDPs can derive both from the insoluble fibrin clot and from soluble fibrin polymers [1]. The

breakdown of both fibrinogen and fibrin by plasmin is counterbalanced by multiple enzymatic modulators (thrombin-activatable fibrinolysis inhibitor, α -2-antiplasmin, and α -2-macroglobulin), whose functions are to limit fibrinolysis at the injury site [3].

The term "fragment D-dimer" was initially used to describe the final plasmin digestion products (resistant to further plasmin breakdown) of the factor XIIIa-cross-linked fibrin clot, which is the fragment DD/E [1,10]. However, the actual D-dimer antigen that can be detected by current immunoassays is not necessarily the DD/E complex [1]. In fact, the term D-dimer comprises a broad mixture of degradation products of cross-linked fibrin with molecular weights likely to range from 190 to over 10,000 kDa [9]. D-dimer fragments are mainly eliminated by renal clearance and reticuloendothelial system catabolism. The plasma half-life of D-dimer is 8 h [11], which is much longer than that of the thrombin-antithrombin complex (10–15 min) or prothrombin fragments 1 + 2 (F 1 + 2; 90 min). In normal physiological conditions, around 2–3% of fibrinogen is converted into fibrin, which instantly enters the fibrinolytic pathway. Thus, D-dimer is measurable in small amounts in healthy subjects and tends to increase with aging [4,12].

Fibrinogenolysis and fibrinolysis have to be considered differently; fibrinogenolysis leads to generation of fragments X, D, Y, and E (fibrinogen degradation products), whilst fibrinolysis leads to more complex fragments called X-oligomers that contain the D and E

fragments in various sequences [1,4]. Only fibrin polymers that have undergone factor XIII-mediated cross-linking will produce fragments containing covalent bonds between two adjacent D domains (i.e. D-dimer). Therefore, D-dimer (DD/E fragment) should be considered as a specific FDP that is recognized by most of the reagent antibodies used in the laboratory assessment of thrombosis [3]. Commercial D-dimer immunoassays are currently predicted to detect an epitope on the degradation products of factor XIIIa-cross-linked fibrin. Monoclonal antibodies employed should not detect fibrinogen degradation products (e.g. X, Y, D, E) on non-cross-linked FDP, but they predictably do detect only epitopes present in the factor XIIIa-cross-linked fragment D domain of fibrin [1,4]. Gaffney et al. found that antigenic determinants of D-dimer assays are a part of the polypeptides in the D-domain, which are conformationally reactive only after factor XIIIa and plasmin have modified the protein [10]. Notably, each monoclonal antibody has its own specificity toward FDP [13].

Preanalytical variables

According to ISO 15189:2012 Medical laboratories—Requirements for quality and competence (ISO, International Organization for Standardization), the ISO standard that is used for laboratory accreditation, the preanalytical phase is defined as “processes that start, in chronological order, from the clinician’s request and include the examination request, preparation and identification of the patient, collection of the primary sample(s), and the transport to and within the laboratory, and end when the analytical examination begins” [14].

Preanalytical phase in hemostasis laboratories

It is now clearly established that most errors encountered in the hemostasis laboratory are related to preanalytical activities [15–20]. Preanalytical errors have a frequency of 60–70%, a frequency that is much higher than that occurring in the analytical (e.g. 10–15%) and postanalytical phase (e.g. 15–20%) [20]. Preanalytical errors are related mainly to intensive manual activities [20,21].

During a 2-year follow-up period, Salvagno et al. studied all preanalytical errors detected in the local coagulation laboratory, and found that the most frequent preanalytical problems were related to samples not received in the laboratory (49.3%), hemolysis (19.5%), clotting (14.2%), and inadequate sample volume (13.7%) [22]. They also showed that preanalytical problems could be identified in up to 5.5% of all coagulation samples [22]. Another study performed by Grecu

et al. during a 1-year period revealed that clotted samples were the major source of preanalytical errors in their hematology laboratory (43.2%), whilst samples with inappropriate volume were less frequent [23]. Recently, Dikmen et al. performed a similar study during a 1-year follow-up and observed that clotted samples (35%) and inadequate volume (13%) were the major causes of sample rejection [24]. They found that the specimen rejection rate for coagulation tests was higher (13.3%) than all other tests (e.g. 3.2% for biochemistry tests, 9.8% for blood gases and 9.8% for urinalysis) [24].

The source of preanalytical variables impairing sample quality can be divided into three main categories: (i) sample collection (e.g. needle size, collection tubes), (ii) sample delivery to the laboratory (e.g. pneumatic tube system (PTS), temperature), and (iii) sample processing (e.g. centrifugation, hemolyzed samples). The storage and stability of samples, as well as the freeze-thaw effect, are also part of the preanalytical process [25].

Importantly, each preanalytical step is vulnerable to errors. Therefore, following preanalytical recommendations at each step of the process is essential for preserving sample integrity [25]. The impact of each of these preanalytical variables on D-dimer measurement in current hemostasis laboratories will be discussed.

Through a review of the literature of the effects of preanalytical variables on D-dimer, it appeared that statistical analysis was frequently used to assess a significant preanalytical change. The absence of significant bias is considered a robust indicator that a pre-analytical variable has not significantly impacted a particular analyte. However, a significant difference does not inherently correspond to clinical significance. Therefore, a clinical criterion is also recommended to accurately assess the effect of preanalytical variables on a particular analyte. Regarding D-dimer, a clinically acceptable cutoff of 10% has frequently been used to assess a significant preanalytical change (especially to assess the effect of interfering substances and for performing stability studies). Unfortunately, this 10% cutoff has not been derived from biological variation data (i.e. reference change values [26]) because such data are still lacking in the hemostasis field. Further validation of this empirical 10% cutoff and establishment of biological variation studies are hence needed [27].

Preanalytical variables affecting the D-dimer measurement

Sample collection

Butterfly devices and needle bore size

It is currently recommended to use straight needles with a diameter ranging from 19 to 22 gauge (G) [21].

Blood samples should be obtained with relatively non-traumatic venipuncture. Excess vein manipulation with the needle should be avoided to limit the risk of developing clots [25], which have been found to potentially impact D-dimer concentration [28]. The use of butterfly devices (small needle attached to flexible plastic wings and connected with a flexing tube) is usually discouraged because the passage of blood along the tubing might be a cause of enhanced shear stress, thus potentially triggering hemostasis activation, erythrocyte injury (i.e. hemolysis); they are also more expensive [21,29,30].

Interestingly, Lippi et al. found that the use of butterfly devices (21 G, 300 mm-long tubing) did not introduce a significant bias in D-dimer values measured with the Mini Vidas® Immunoanalyzer (bioMérieux, Marcy l'Etoile, France) compared to conventional straight needles (21 G) [31]. This observation implies that the potential activation of the hemostatic system occurring within the tubing of the butterfly device, regardless of its composition, may be negligible for the D-dimer measurement.

Lippi et al. also found a clinically negligible bias using butterfly devices with different needle bore sizes (21, 23, or 25 G), when D-dimer was again measured with the Mini Vidas® Immunoanalyzer [32]. Taken together, the use of butterfly devices, even with small size needles, seems to represent a viable option to standard straight needles [31,32]. However, straight needle venipunctures are preferable, especially when other coagulation tests are requested [33]. The use of butterfly devices may, therefore, be considered in particular individuals (geriatrics, oncology, pediatrics, or in emergency settings) [21,33]. When blood is collected using a winged collection system (or IV catheters), a discard tube is mandatory for removing air contained within the tubing, which may be associated with the collection of an inadequate volume of blood [34].

Tube material

Non-activating material (silicone-coated glass or polypropylene plastic) is preferred for hemostasis testing because this will prevent initiation of clotting due to spurious hemostasis activation in the blood collection tube [25].

Several studies were performed for comparing the use of different collection tube materials. Leroy-Matheron et al. found no significant differences in measured D-dimer values (Asserachrom® ELISA assay) using glass or polyethylene terephthalate (PET) plastic collection tubes (Vacutainer®, Becton Dickinson (BD)) [35]. Gosselin et al. also failed to find significant differences for D-dimer values (Advanced D-Dimer) when using plastic (Vacutainer®,

BD and Vacuette®, Greiner Bio-one (polypropylene)) and glass citrated collection tubes (3.2%) [36]. More recently, Yavas et al. compared three plastic citrate (3.2%) tubes (Vacutainer®, BD; Vacutainer® Plus Plastic, BD; and Vacuette®, Greiner Bio-one) with a standard glass tube (Vacutainer®, BD), and also failed to find significant differences [37]. Therefore, the quality of D-dimer measurement seems to be not substantially dependent on the type of collection tube.

Anticoagulant sample (or tube additive)

The use of collection tubes containing 3.2% (105–109 mmol/L) buffered sodium citrate anticoagulant is now recommended by the Clinical Laboratory Standards Institute (CLSI) and by the World Health Organization (WHO) for the vast majority of hemostasis tests [25,34]. A blood to anticoagulant ratio of 9/1 is recommended because the sodium citrate anticoagulant can be used only in a liquid form [25]. Failure to correctly fill the citrate anticoagulant tube will typically prolong clotting times (PT, APTT, and TT), and may also lead to underestimating D-dimer and fibrinogen [28]. Therefore, the sample should be checked for quality (e.g. clotting, underfilling, or overfilling) [4]. Blood should also flow freely into the tube and be promptly mixed within 30 s after phlebotomy (i.e. from 3 to 6 complete inversions), thus ensuring the complete distribution of anticoagulant activity [25]. Serum and heparinized/EDTA plasma samples cannot be accepted because these additives dramatically interfere with clot-based assays (e.g. PT, APTT, FV, FVIII) [25,28].

However, several D-dimer assays (including POC tests) allow the use of citrated, heparinized or EDTA plasma (Pathfast®, Tina-quant® (Roche, Switzerland)), AQT-90® (Radiometer, Denmark), Simplify® (Agen Inc, Australia)), whilst others recommend only citrate anticoagulated blood (e.g. Vidas®, BCS®, STA-Liatest® (Diagnostica Stago, Asnières, France), Immulite® (Siemens, USA)) [38]. Schutgens et al. observed a higher mean D-dimer concentration (Tina-quant® D-dimer) in heparin collection tubes (2510 µg/L), compared to citrate (2060 µg/L) albeit the difference did not achieve statistical significance [39]. The higher D-dimer concentration in heparin plasma is probably attributable to the introduction of a dilution factor in citrate-containing blood tubes. Although Vukovich et al. reported that heparinized and citrated samples yielded identical D-dimer values (D-dimer Gold®, Agen Biochemical) samples, they recommended using a correction factor of 0.84 to take into account the dilution of the citrate anticoagulant, a process that may be confusing in routine practice [40]. Local evaluation of other matrices should be validated

before they are used in clinical practice, especially when the manufacturer recommends using only citrated samples. Lippi et al. found only a modest bias of D-dimer values (Immulate[®] 2000 D-dimer) between lithium-heparin and sodium citrate blood tubes [41].

The main advantage of using other anticoagulants (e.g. heparin) is the possibility to measure other analytes (e.g. electrolytes, cardiac troponins). However, the dilution factor may be a confounding factor and may require setting specific algorithms for correcting data. Therefore, other anticoagulants should be locally validated or their use should be reserved for special circumstances; thus buffered sodium citrate (3.2%, or 105–109 mmol/L) remains the recommended sample matrix [33].

Notably, D-dimer was first measured in serum samples with preliminary removal of fibrinogen to avoid cross-reactivity with polyclonal antibodies [4]. However, false-positive results were frequently encountered in patients under anticoagulant treatment, whilst false negative values were also seen when FDPs were entrapped in the clot. The loss of fibrinopeptide A has also been associated with degradation of non-cross-linked FDP [4].

Tourniquet use

Tourniquets are conventionally used to temporarily obstruct the vein flow and thereby assist the phlebotomist to identify vein access. It is commonly recommended that the tourniquet should be removed as soon as the needle is in the vein or when the first tube starts to fill and that it should never remain in place for more than 1–2 min [21,25]. The most important drawbacks of prolonged tourniquet placement include hemoconcentration and clot formation, which may jeopardize the quality of coagulation testing [21,25,32]. More specifically, Lippi et al. observed that D-dimer values (Mini Vidas[®] Immunoanalyzer) were significantly increased, by 13.4%, when measured in samples collected after 3 min of venous stasis [32]. Lower differences were also observed after a 1-min stasis (mean increase of 7.9%) [32].

Sample delivery to the laboratory

Samples should be delivered to the laboratory at ambient temperature (15–22°C), in the shortest possible time (usually <1 h) after collection [34]. The stability of D-dimer will be discussed below. Blood samples have to be delivered in a vertical position and the cap should not be removed [21,25]. The transport of blood tubes in a vertical rather than a horizontal position limits the

extent of *in vitro* microparticle generation [21,42]. Although the presence of blood collection facilities close to the laboratory is the best option to limit pre-analytical variability, many hospitals deliver blood samples with a PTS [21]. PTS has the great advantage of reducing the turnaround time (TAT), but these systems need to be validated before use, because excessive acceleration/deceleration, radial gravity forces, vibration, and changes in air pressure may trigger platelet activation and hemolysis [21,43,44].

Schutgens et al. failed to find significant difference between D-dimer values (Tina-quant[®] D-dimer) measured in samples delivered by PTS (≈100 m long) compared to those carried by hand [39]. Wallin et al. also found no substantial difference in D-dimer values (MediRox[®] assay, Nyköping, Sweden) delivered to the laboratory by PTS (≈500 m long) compared to those directly collected in the laboratory [44]. Unlike these studies, Le Quellec et al. observed significant difference of D-dimer values (DDi-HS 500[®], HemosIL[®]) in samples conveyed by PTS (≈2000 m long) compared to those transported with a motor vehicle (mean difference of 7.4%) [43]. Nevertheless, a disagreement at the 500 µg/L D-dimer cutoff could be observed only in 1 out of the 39 samples (510 µg/L with PTS and 490 µg/L by motor vehicle) and the overall difference was still within the intra-assay imprecision of the D-dimer assay. Although the data in the literature seemingly is reassuring about the low impact of PTS on D-dimer test results, it is advisable that each laboratory assesses its local PTS because systems are heterogeneous in terms of length, internal diameter, maximal acceleration force, and speed [43].

Sampling processing

Prior to being analyzed, samples should be carefully checked for the presence of inappropriate additive, identification errors, insufficient volume (9/1 ratio not respected), and even for the presence of clots. Unsuitable samples should then be rejected and test results suppressed to preserve patient safety [21].

Centrifugation

Except for point of care (POC) whole blood D-dimer analyzers, samples are centrifuged at room temperature at 1500×g for at least 15 min, so that plasma can be reliably separated from the cellular components [34]. According to Bernard et al., however, no difference could be observed in D-dimer values (Vidas[®] assay) when samples were centrifuged at 4500×g for 2 min rather than at 1500×g for 15–20 min [45]. The use of a shorter centrifugation time will be indeed useful in

emergency settings to reduce the TAT. The samples can be kept at room temperature for up to 8 h and at 4 °C for up to 24 h (unspun citrated blood samples) before measurement of D-dimer [4,46,47]. Elf et al. centrifuged blood samples at 4 °C and measured D-dimer with two quantitative immunoassays (Innovance® and AxSYM® (Abbott, USA)) [48], although there was no evidence supporting cold centrifugation for D-dimer testing.

Interfering substances

The types of interferences occurring in the preanalytical phase of coagulation testing are typically classified into paraproteinemia, icterus, lipemia, and, the most studied, hemolysis [20].

In vitro hemolysis still represents one of the most frequent causes of preanalytical problems in clinical laboratories, with a prevalence ranging between 30–70% of all unsuitable specimens [22,49–52]. Causes of hemolyzed samples could result from patient's characteristics or diseases (e.g. hemolytic anemia, metabolic disorders, infectious agents, hemoglobin-based blood substitute), phlebotomy (e.g. needle G, tourniquet time, traumatic draw, no mixing or vigorous mixing of blood tubes), sample transport (e.g. PTS, time, emergency or intensive care units), processing (e.g. delay before centrifugation, specimen re-spun, time/temperature of transport), and storage (e.g. temperature and duration) [50,51]. The interference will also depend on the analytical method used (photometric, clotting, immunoassay). The CLSI guideline H21-A5, collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays, advises against analyzing samples with visible hemolysis because of possible bias in coagulation tests from release of pro-coagulant factors from injured cells [34].

Lippi et al. showed that a significant increase in D-dimer level (Vidas® assay) could be observed in samples containing a final whole blood lysate concentration of at least 2.7% obtained with a freeze-thaw cycle (–70 °C) [53]. However, a clinically significant variation (10%) was observed only above cell-free hemoglobin concentrations of 13.6 g/L (or final lysate concentration of ≈6.4%) [51,53]. Another study by the same authors evaluated the impact of mechanical hemolysis (syringe equipped with a fine needle (30 G, 0.3 × 8 mm)) on two different D-dimer assays (HemosIL AcuStar® (Werfen, USA) and HemosIL HS® (Werfen, USA)) [49]. The D-dimer concentration measured with the HemosIL HS® reached a clinically significant decrease (–5%) from a hemolysis index between 5.5–7.0 g/L of cell-free hemoglobin. Greater hemolysis (i.e. 11.5–15.0 g/L of cell-free

hemoglobin) was necessary to achieve the same level of D-dimer decrease (–5%) with the AcuStar D-dimer assay [49]. However, the bias never exceeded the clinically acceptable cutoff of 10%. D'Angelo et al. more recently reported a significant difference of D-dimer values with mild hemolysis (cell-free hemoglobin of 0.5 g/L) obtained by either heat-shock at –80 °C (+6.2%) or by mechanical injury (rotating blade homogenizer for 30 s) (+14.9%) when measured on a BCSxp® assay (Siemens Healthcare) [54].

The comparison of studies assessing hemolysis interference is challenging. First, different D-dimer assays with different reagent formulations were used. Secondly, different techniques to obtain hemolyzed specimens were used (freezing and thawing whole anticoagulated blood, mechanical lysis of whole anticoagulant blood by aspiration through a fine needle or by a rotating blade homogenizer). While spiking samples with hemolysate or pure hemoglobin solution is considered unsuitable (lysis of leukocytes and platelets is taken into account), mechanical hemolysis obtained with a syringe equipped with a fine needle is preferred (i.e. it more closely reproduces the breakdown of cells due to traumatic blood collection) to the freeze-thaw method, which also suffers from a lack of standardization across studies (e.g. different temperatures (–70 °C [53] or –80 °C [54]) and duration of freezing [51]. Thirdly, the use of a clinical cutoff derived from biologic variability studies should be preferred to the analytical bias assessed [51]. Fourth, the selection of blood samples may produce different results and conclusions; data generated in healthy volunteers may not necessarily be similar to that obtained in intensive care unit patients, or in those undergoing multiple therapies (namely, antiplatelet or anticoagulant drugs) [51].

The rejection of all hemolyzed samples is recommended by the CLSI. The majority of hemolyzed samples (±95%) from clinical laboratories are only mildly hemolyzed (cell-free hemoglobin 0.3–0.6 g/L) [49,50,55]. Thus, when the concentration of cell-free hemoglobin is below a non-interference limit (i.e. <3 g/L), D-dimer values may still be reliable and may be safely reported to the clinicians [49]. Notably, avoiding rejection of all hemolyzed specimens will reduce additional blood sampling and may shorten the clinical decision-making process, thereby favorably impacting staff workload, patient comfort, and costs [49,51].

The effect of lipemia, icterus, and paraproteinemia has been less widely discussed in the literature [51,56,57]. In an early report, Pittet et al. failed to find any impact of lipemia (triglycerides concentration of 9.7 g/L) and icterus (bilirubin concentration of 485 mg/L,

829.5 $\mu\text{mol/L}$) on the Vidas D-dimer assay [58]. More recently, La'ulu et al. found similar findings on the AxSYM D-dimer assay, with no clinical significant difference ($<10\%$ of deviation from the baseline) following addition of bilirubin (final concentration of 292 mg/L, 499.4 $\mu\text{mol/L}$) and triglycerides (final concentration of 41.6 g/L) [59].

The Innovance D-dimer assay on the Sysmex CA-7000 also appeared to be unaffected by bilirubin (both free and conjugated bilirubin up to 1000 mg/L, 1710.4 $\mu\text{mol/L}$) [60]. However, spurious D-dimer results were observed on the STA analyzer after adding high concentrations of free bilirubin (from 800 mg/L, 1368.3 $\mu\text{mol/L}$) [60]. According to Chen et al., D-dimer, as measured with the Sysmex CS-5100 analyzer, was also free from triglycerides and total bilirubin interferences [61]. Moreover, the HIL check was able to detect lipemic samples (triglycerides concentrations >4 g/L) on the CS-5100 D-dimer assay in three out of four samples [62]. The undetected sample had a very high concentration of triglycerides (9.43 g/L).

Regarding the paraproteinemia interference, Mugler et al. reported a case of a patient with Castleman disease associated with an elevated D-dimer level [57]. During investigation of this finding, no D-dimer band could be noticed on Western blot compared to D-dimer measured with the STAR coagulation instrument. The patient's Castleman disease-associated monoclonal gammopathy was identified to be the source of the false positive D-dimer result. Paraproteinemia interference was also hypothesized in another report published by Roller et al. and it could not be excluded from the reports of Wu et al. and Huang et al. [63–65].

Spurious D-dimer concentrations due to heterophilic antibodies have also been reported in the literature [56,66,67]. Potential interferences in D-dimer testing should always be suspected whenever clinical discrepancies arise. Heterophilic blocking reagents could be used to overcome the problem [56,57]. The comparison with another method to point out interferences represents also a good choice to search for interference [57,68].

Stability, storage, and freeze-thaw effects

It is currently recommended that samples should be maintained at room temperature (15–25 °C) for no more than 4 h before testing [21]. However, and as discussed below, several studies have shown that D-dimer may be stable for a longer period. This is of particular concern because clinicians frequently ask to add the

D-dimer analysis to samples already delivered from the emergency room.

Detailed information on D-dimer stability is shown in Table 1. With the exception of the publication of Schutgens et al., citrated specimens were used to assess the stability of this analyte [39]. The D-dimer stability was studied extensively in both plasma and whole blood. In brief, D-dimer is stable for at least 24 h (in plasma/whole blood) at room temperature (RT) or at 2–8 °C with many different D-dimer immunoassays. Only Böhm-Weigert et al. and Toulon et al. suggested, respectively a 6- and 8-h D-dimer plasma stability, but this was because they did not extend their analysis for longer periods [47,69]. Frozen samples may also be used when longer storage periods (i.e. months or years) are needed, especially for research purposes [69–72].

Several criteria are available to assess stability of analytes [73]. As discussed above, statistical analysis is frequently used, and the absence of significant bias is considered a robust indicator of stability. However, a significant difference does not inherently correspond to clinical significance [74]. Therefore, a clinical criterion is also recommended to accurately assess stability of a particular analyte. It has also been demonstrated that the use of different stability criteria may strongly impact the result of stability studies [73]. Therefore, the use of multiple criteria is encouraged. All studies published so far used multiple stability criteria except that performed by Betsou et al. [70]. The review of D-dimer stability studies showed that a clinical cutoff of 10% was most frequently used [69,71,72,75–78]. However, this 10% cutoff has not been derived from biological variation data. Further validation of this empirical 10% cutoff is needed, even though no discordance could be found in the statistical approach available in the literature. The analysis of discordance based on the cutoff used to rule out venous thromboembolism (VTE), as evaluated by Caliezi et al., Böhm-Weigert et al., and Toulon et al. [47,69,79], is also a criterion that seems reasonable.

The studies of Schutgens, Zürcher, and Gosselin also showed that the freeze-thaw procedure did not significantly affect D-dimer concentration measured with three different D-dimer assays [39,75,80]. Böhm-Weigert et al. evaluated the impact of four freeze-thaw cycles (-60 °C or less) and failed to find clinically significant changes of plasma D-dimer values (mean deviation always $<10\%$) [69].

In conclusion, while D-dimer is mainly used for emergency diagnoses and when rapid measurement is needed, stability information is required when samples are stored in central laboratories and when they are

Table 1. Stability of D-dimer.

Stability	Conditions	Anticoagulant	Plasma/whole blood	D-dimer assay	Subjects	Stability criteria	Reference
24 h	RT	Heparin	Plasma	Tina-quant [®] (Roche)	17 patients	Student t-test and regression equation	[39]
24 h	RT	Citrate	Plasma	Tina-quant [®] (Roche)	15 patients	Student t-test and regression equation	[39]
6 h	RT	Citrate	Plasma	Innovance [®] (Siemens)	40 patients	10% deviation from baseline, regression equation and discordance at the cutoff level of 0.5 mg/L FEU	[69]
24 h	RT	Citrate	Plasma	Innovance [®] (Siemens) ^a	80 patients	10% deviation from baseline, analysis of variance, regression equation and Pearson correlation coefficient	[78]
24 h	RT	Citrate	Whole blood	Vidas [®] (bioMérieux)	117 patients	Spearman correlation coefficient, regression equation and discordance at the cutoff level of 500 ng/mL FEU	[79]
24 h	RT	Citrate	Whole blood	Innovance [®] (Siemens)	44 patients	10–20% deviation from baseline, Student t-test, and regression equation	[76]
24 h	RT	Citrate	Whole blood	ACL-TOP [®] (Werfen)	26 patients	Wilcoxon's paired t-test, regression equation, and bias plot	[46]
52 h	RT	Citrate	Whole blood	Asserachrom [®] (Stago)	59 patients	Analysis of variance, 10% deviation from baseline	[75]
8 h	RT	Citrate	Whole blood	ACL-TOP [®] (Werfen)	144 patients	Analysis of variance, Student t-test or Wilcoxon signed rank test, Bland-Altman plot and discordance at the cutoff level of 0.5 ng/mL FEU	[47]
24 h	2–8 °C	Citrate	Plasma	Innovance [®] (Siemens)	40 patients	10% deviation from baseline, regression equation, and discordance at the cutoff level of 0.5 mg/L FEU	[69]
24 h	4 °C	Citrate	Plasma	Innovance [®] (Siemens) ^a	80 patients	10% deviation from baseline, analysis of variance, regression equation and Pearson correlation coefficient	[78]
24 h	4 °C	Citrate	Plasma	Vidas [®] (bioMérieux)	20 patients	Wilcoxon's paired t-test, 10% deviation from baseline	[77]
24 h	4 °C	Citrate	Whole blood	ACL-TOP [®] (Werfen)	26 patients	Wilcoxon's paired t-test, regression equation, and bias plot	[46]
24 months	–24 and –75 °C	Citrate	Plasma	STA-Liatest [®] (Stago)	Plasma pool (6 patients)	Statistical change ^b , 5–10% deviation from baseline	[72]
2 weeks	–20 °C	Citrate	Plasma	STA-Liatest [®] (Stago)	23 hV and 18 patients	Paired t-test, 10% deviation from baseline	[71]
36 months	–60 °C (or less)	Citrate	Plasma	Innovance [®] (Siemens)	40 patients	10% deviation from baseline, regression equation, and discordance at the cutoff level of 0.5 mg/L FEU	[69]
9 years	–80 °C	Citrate	Plasma	STA-Liatest [®] (Stago)	60 patients	Wilcoxon's paired t-test	[70]

RT: room temperature; FEU: fibrinogen equivalent units; hV: healthy volunteers.

^ausing the Sysmex[®] CA 7000 platform.^bspecific test not mentioned.

Table 2. Specific pre-analytical data regarding D-dimer testing.

Pre-analytical variables	General recommendations in hemostasis laboratories	Specific data regarding D-dimer
Sample collection		
Needle bore size	19–22 G	23–25 G also tolerated
Butterfly devices	Discouraged	Tolerated
Tube material	Non-activating material (silicone-coated glass or polypropylene plastic)	Glass or plastic
Anticoagulant sample	Sodium citrate 3.2% (105–109 mmol/L)	Sodium citrate or heparin*
Tourniquet use	Removed as soon as the needle is in the vein (max 1–2 min)	Longer tourniquet use (i.e. 3 min) not tolerated
Sample delivery to the laboratory	At RT (15–22 °C), in vertical position, usually <1 h	PTS tolerated
Sample processing		
Centrifugation	At RT, 1500×g for at least 15 min	Faster protocol allowed (at RT, 4500×g for 2 min)
Interfering substances	Do not analyze samples with hemolysis	Cell-free hemoglobin i.e. <3 g/L tolerated
Stability, storage and F/T effects	At RT (15–22 °C), no more than 4 h	At least 24 h at RT or at 2–8 °C or years at –60 to –80 °C
		No impact of F/T procedure

RT: gauge; RT: room temperature; PTS: pneumatic tube system; F/T: freezing/thawing.

*correction factor needed (dilution).

used in a research environment [33]. Table 2 summarizes specific pre-analytical data regarding D-dimer testing.

Analytical variables

D-dimer assays

The first generation of D-dimer assays, which were developed in the 1970s, was capable to detect both fibrinogen and FDP by means of polyclonal antibodies [2]. These immunoassays could be performed only in serum to avoid the cross-reactivity of fibrinogen present in high concentrations in plasma. However, false-positive results were encountered in patients taking anticoagulants, and false negative results were encountered in samples that developed clots and when degradation products absorbed into the clot.

The loss of fibrinopeptide A during serum preparation impairs detection of non-cross-linked FDP [4]. The measurement of FDP was initially performed using (among others) methods based on latex fixation and agglutination, hemagglutinin inhibition, staphylococcal clumping, immunoelectrophoretic, and immunodiffusion [3].

A significant gain of specificity and sensitivity was reached with the introduction of monoclonal antibodies targeting D-domains in the early 1980s [3]. These immunoassays allowed the specific targeting of D-dimer epitopes, which were lacking in both FDP and non-cross-linked fragments of fibrin [4,81]. The use of plasma was, therefore, possible because of minimal fibrinogen (or FDP) cross-reactivity. Since the first monoclonal antibody developed by Rylat et al. in 1983 (3B6), more than 20 monoclonal antibodies have been produced and used in laboratory practice [9,81]. First generation methods were represented mainly by qualitative latex agglutination immunoassays using anti-body-coated latex microparticles and needing visual

inspection of agglutination [1,3,4,82]. The development of automated second-generation latex agglutination immunoassays (or the latex-enhanced immunoturbidimetric assay) allowed quantitative measurement of D-dimer through recording aggregation rate in response to D-dimer [1,2]. Once D-dimer was added, latex microparticles agglutinated, thus preventing the light to pass through the solution. The increase of light absorbance measured with a photometer was directly proportional to the D-dimer concentration [3]. Latex-enhanced immunoturbidimetric assays are rapid and achieve comparable sensitivity to conventional enzyme-linked immunosorbent assays (ELISA) [3].

Microplate ELISA assays were employed before latex agglutination assays, mainly for research purposes [1]. In these assays, a capture antibody binds D-dimer antigen on a plate. After incubation, a labeled antibody is added to the well and binds the immobilized D-dimer antigen (sandwich assay) to promote a colorimetric reaction [1,2]. Although microplate ELISAs are very sensitive to D-dimer and have long been considered the reference method [1–3], they are manual, require technical skills, are time-consuming (i.e. around 2–4 h), and are plagued by a high degree of analytical imprecision [2,83].

In the mid-1990s, bioMérieux developed an ELISA-based assay with fluorescence end-point detection (enzyme-linked immunofluorescence assay (ELFA)) [3,58]. This assay displayed sensitivity and specificity similar to those of microplate ELISAs, with the great advantage of being automated, thus generating more precise results in a shorter time (30 min) [3,58]. The Vidas® assay is still considered the reference commercial quantitative immunoassay [84–86] and is the most clinically validated D-dimer measurement technique [38,48,85,87–89].

More recently, chemiluminescent enzyme immuno-metric assays that display similar sensitivity to ELISAs and latex-enhanced immunoturbidimetric assays have

been developed [3,41,84]. Magnetic particles coated with monoclonal antibodies specific to D-dimer are used; incubation of anti-D-dimer antibody labeled with isoluminol generates a chemiluminescent reaction that is directly proportional to the D-dimer concentration [3].

Utilization of POC D-dimer assays is especially attractive for general practitioners (GP) because quantitative D-dimer assays are not always readily available in rural areas, during weekends, or at home visits [90–92]. In such settings, the GP had to refer the patient to a central laboratory facility [91]. The main interest in POC D-dimer assays is the possibility to rapidly screen patients for thromboembolic disease, thus leading to a decrease in overcrowding in urgent care facilities [3,90,93,94]. The majority of POC tests use whole blood, are homogenous, use monoclonal antibodies and have a short TAT [3,90]. In hemagglutination assays (e.g. SimplyRED®), bivalent antibodies are specific to red blood cells and D-dimer. In the presence of D-dimer in the blood, red cell agglutination occurs and provides a qualitative result [3,4]. Detection based on immuno-chromatography (e.g. Clearview Simplify®), fluorescence (e.g. Vidas® mini analyzer, Triage®, Stratus CS®), and chemiluminescence (e.g. Pathfast®) also exist. A few semi-quantitative assays remain in use (e.g. Dade Dimertest®, NyoCard®) [3,95]. The recent review of Riley et al. compiled the different POC D-dimer assays available on the market [3]. Table 3 summarizes the main characteristics of D-dimer assays.

Inter-laboratory variation

Variability across assays

The fibrin assay comparison trial (FACT) study, published in 2001, evaluated 23 quantitative D-dimer assays (15 latex-enhanced immunoassays, 6 ELISAs, and 2 membrane-based immunoassays) [96]. The authors found that the mean values obtained in 39 samples varied from 630 µg/L to 13,350 µg/L (≈21-fold), with two assays displaying significant cross-reactivity toward fibrinogen degradation products. It was also found that ELISAs and latex-enhanced immunoassays were more reactive to low-molecular-weight cross-linked fibrin and high molecular weight fibrin, respectively [96]. Accordingly, the study of Meijer et al., based on 357 laboratories using the 7 frequently-used D-dimer immunoassays, found that some assays gave 20 times higher D-dimer concentrations compared to others [97]. Another study from 423 laboratories comparing D-dimer values also showed high variability in results close to the cutoff used for VTE exclusion [98]. In 2014, in a survey involving 3800 laboratories, the Coagulation

Table 3. Characteristics of D-dimer assays [2,3,41,84,91,165,174].

Type	ELISA	ELFA	Unenhanced latex agglutination assay	CLIA	Latex-enhanced immunoturbidimetric assay	POC assay
TAT	Quantitative	Quantitative	Qualitative/semi-quantitative	Quantitative	Quantitative	Qualitative/quantitative
Pros	2–4 h Considered as the gold standard, sensitivity, observer-independent	35–40 min Considered as reference method, most validated method, sensitivity, automation, wider linear range, automated, observer-independent	Rapid Rapid, inexpensive	25–40 min Sensitivity, rapid, automated, observer-independent	15 min Sensitivity, automated, rapid, observer-independent	2–20 min Readily available, fast, higher specificity, whole blood
Cons	Highly manual, technical skills, time-consuming, not optimal linear range, moderate specificity	Moderate specificity	Moderate sensitivity, manual, observer dependent	Lack clinical validation, moderate specificity	Moderate specificity	Sensitivity, not all FDA cleared, observer dependent, manual
Examples	Asserachrome® (Stago), Enzygnost® (Dade Behring)	Vidas® (bioMérieux)	Dimertest latex® (IL), Fibrinosticon® (bioMérieux), Dade Dimertest® (Siemens)	AcuStar® (Werfen), Immulite® (Siemens)	Tina-quant® (Roche), STA-Liatest® (Stago), HemosIL HS® (Werfen) Innovance® (Dade-Behring)	Simplified® (Agen), Clearview® (Agen), Simplify® (Agen)

ELISA: Enzyme-linked immunosorbent assay; ELFA: Enzyme-linked immunofluorescence assays; CLIA: Chemiluminescent enzyme immunoassay; POC: Point of care.

Resource Committee of the College of American Pathologists (CAP) found that the inter-method coefficient of variation (CV) was as high as 42% [99]. Based on four UKNEQAS reports (April 2017, July 2017, September 2017, and January 2018), the inter-method CVs (after the exclusion of outliers) were reported to be 33.4%, 38.2%, 41.8%, 30.9% and 18%, 19.3%, 17.7% and 19% for non-FEU and FEU units, respectively; these four surveys included between 736 and 752 laboratories.

Calibrators

The controlled lysis of fibrin clots is used mainly to obtain calibration materials [9]. Hence, manufacturers have to ensure that lysis is reproducible in order to maintain the same variety of size degradation products, because the sensitivity of the assay may change according to the relative amounts of high molecular weight fibrinogen (HMWF) or low molecular weight fibrinogen (LMWF) [9,96]. D-dimer assays may be calibrated with purified fibrin fragment D-dimer equivalents or according to the fibrinogen amount used to prepare the calibrator [9]. D-dimer units thus may vary according to the type of calibrator used, i.e. D-dimer units (DDU) or fibrinogen equivalent units (FEU) [9,96,100].

Standardization versus harmonization

The heterogeneity of fragments derived from plasmin digestion of cross-linked fibrin (from LMWF to HMWF) [9], the use of monoclonal antibodies with different specificities towards D-dimer epitopes [1,96,101,102], the lack of international certified internal controls or calibrators [103] and the use of different units and clinical cutoffs are the leading sources of the large inter-laboratory variability, thus making D-dimer assay standardization a challenging, if not impossible, target [9,13,99,101,102,104–108]. Because the requirements for D-dimer assay standardization cannot be met, less stringent harmonization procedures, based on the use of mathematical models, have been proposed to make results obtained with different assays more comparable [96,97,100,102].

The first attempt to achieve D-dimer assays harmonization was published in 1997 by Nieuwenhuizen et al. [102]. Different pools of patients with various diseases were prepared and tested with 5 different D-dimer assays (one microlatex and four microplate ELISA), and a "pool consensus value" (corresponding to the mean D-dimer result of each assay) was assigned to each pool. Results of each pool obtained from the individual assays were then divided by the corresponding "pool consensus value" and averaged. Samples used to

prepare pools were also measured to verify that this model could be broadly implemented. The squared regression coefficient values ranged from 0.7 to 0.92, which led the authors to conclude that the use of a conversion factor may be a feasible approach to harmonized data generated with different assays [102].

The second attempt was published by Dempfle et al. in 2001. Overall, 39 individual samples with 23 D-dimer assays (including microlatex-enhanced, membrane-based and ELISA assays) were tested [96]. For each assay, a conversion factor was calculated by using the median value for each sample measured with all assays, and for each assay, the median value obtained with all samples. The multiplication of individual sample assay values by the corresponding conversion factor was found to be effective to improve the correlations among most assays.

The third attempt was conducted by Meijer et al. in 2006 [97]. A plasma pool of 50 patients diluted with normal plasma was used to prepare five different samples that were then distributed in an external quality control survey to 502 participants using seven different D-dimer assays. For each D-dimer assay, the mean results of each sample were plotted against the amount of pool added (assay-specific regression line). The median results of the mean results obtained with each D-dimer assay were also plotted against the amount of pool added (reference regression line). Harmonized results could, therefore, be calculated from the measured result by using the intercept and the slope of both regression lines (assay-specific and reference). The authors showed that a significant decrease in between-assays variation was observed by introducing this model.

The fourth attempt was published by Jennings et al. in 2007 [100]. Three calibrators and two test samples were delivered to more than 500 laboratories that used 9 different D-dimer techniques and that participated in the UK NEQAS external quality survey. Individual laboratory results for calibrators were plotted against the median results obtained with all the D-dimer immunoassays. The individual regression line was used to convert data generated on the two test samples into harmonized results. This approach was effective in improving the between-center agreement after calibration, with significant improvement in inter-laboratory variability for the two samples (from 25.9% to 11.6% and from 22.4% to 7.7% for results reported as FEU; from 55.3% to 21.6% and from 40.8% to 11.6% for results reported as DDU).

Harmonization of D-dimer assays may be particularly useful for diagnosis and monitoring of disseminated intravascular coagulation (DIC), because the high inter-variability between D-dimer assays [96,98] may directly

impact the DIC score and lead to misclassifications. Harmonization of D-dimer results would provide better consistency and improve the calculation of DIC scores [9]. However, the inherent heterogeneity and the challenges of the harmonization efforts make the identification of a unique diagnostic cutoff even more difficult and clinically questionable, because false-negative results have an unfavorable impact on patient management [9].

Nevertheless, the above-mentioned approaches indicate that harmonization is feasible, and it should be considered on a larger scale [103]. Further validation of these models with more assay systems is needed to confirm whether such models may be generalized to all D-dimer immunoassays [9,97,100]. Recently, a call for harmonization was raised by the Fibrinolysis and DIC Standardization Subcommittees of the ISTH [103]. Utilization of a large number of samples from patients with various clinical conditions may be useful to produce a stable freeze-dried reference material containing high concentrations of D-dimer derived from both LMWF and HMWF [103]. A range of D-dimer values could be generated from this reference material, to obtain a reference regression line for D-dimer immunoassays. The use of monoclonal antibodies displaying high affinity to low and middle molecular weight FDP species combined with monoclonal antibodies targeting higher molecular weight forms may also represent a feasible alternative [101]. Continuing discussion among manufacturers, scientists, and clinicians is essential to achieve better harmonization.

Recommendations on the performance of D-dimer assays

Because of the heterogeneity of the different D-dimer immunoassays (i.e. format, antibodies) and the lack of standardization/harmonization, the analytical and clinical performances of the local method should always be evaluated before implementation in management strategies for VTE [109,110].

The cutoff value used to exclude VTE needs to be confirmed locally, using, as suggested by the British Committee for Standards in Haematology guidelines, a minimum of 200 subjects [1,111]. Several manufacturers also recommend local validation of cutoffs [1]. However, this approach may not be practical for all laboratories; thus the cutoff proposed by the manufacturer may be used [4], provided that reliable validation studies have been conducted elsewhere and that no significant batch-to-batch variability has been found [9]. Cut-offs have been clinically validated in prospective studies for a number of assays (e.g. Vidas[®], AxSYM[®],

STA-Liatest[®]). This may be a consideration when choosing a D-dimer assay [1,4,48,85,87,88,107,112]. Otherwise, comparison with validated assays is encouraged. Manufacturers should stay abreast of the recent literature regarding the use of their immunoassays and update their cutoffs when indicated [113].

A CAP survey showed that 488 laboratories out of 1506 in the USA were using cutoff values higher than those recommended in the literature or by the manufacturer [99]. A European survey also highlighted that 24% and 55% of participants used lower or higher cutoffs, respectively, than those recommended [9].

A recent Italian consensus document recommends the use of certified, quantitative immunoassays for testing the plasma of emergency department patients [33]. It was also recommended that D-dimer immunoassays should have an imprecision <10% close to the diagnostic cutoff and a measurement range and linearity of 50–5000 µg/L FEU. Importantly, D-dimer assays should not react with fibrinogen or FDP, and preferably should not react with fibrin and fibrinogen fragments released from proteolysis mediated by various enzymes [103].

Clinicians need to be aware of the performance characteristics of the D-dimer assay they are using because D-dimer comprises a mixture of degradation products of cross-linked fibrin with heterogeneous molecular weights [1,9,105]. Moreover, the laboratory should always use clinically validated cutoff values, because this diagnostic threshold plays a crucial role in clinical decision making [1].

Postanalytical variables

According to the ISO 15189:2012 standard, the postanalytical phase is defined as "processes following the examination including systematic review, formatting and interpretation, authorization for release, reporting and transmission of the results, and storage of samples of the examinations" [14].

Although it is now clearly established that most problems in the hemostasis laboratory emerge from the preanalytical phase (60–70%), [20] the postanalytical phase is still the cause of a large number of diagnostic errors (15–20%) that can jeopardize patient safety [105]. Here we discuss postanalytical issues that are specific to D-dimer measurement.

D-dimer units

Two units of D-dimer measurement are currently used in clinical laboratories: FEU and DDU. On the one hand, FEU compares the mass of D-dimer to that of

fibrinogen, and the calibrators are prepared from plasmin degradation of purified fibrinogen clotted in the presence of factor XIII. On the other hand, DDU determine the mass of the estimated weight of D-dimer, and the calibrators are composed of purified D-dimer [4]. FEU (340 kD) and DDU (195 kD) may be used interchangeably if multiplied by 2 as FEU is roughly two times the mass of one DDU [2–4,9,33]. Another difficulty arises from using as many as 7 different measure units (i.e. ng/mL, mg/L, $\mu\text{g/L}$, $\mu\text{g/mL}$, g/L, g/mL, and mg/dL) for reporting D-dimer results [3]. Furthermore, the same D-dimer assay is currently reported in different units [2,114,115]. Therefore, at least 14 combinations for D-dimer measurement coexist.

As recently studied by Olson and Lippi, the majority of clinical laboratories (59% and 60%, respectively) used FEU for D-dimer measurement, whilst the leading measure unit was mg/L, followed by ng/mL, in both surveys [99,105]. Among the measure units that can be adopted, " $\mu\text{g/L}$ " (or "ng/mL") is probably the unit that best approximates the *Système Internationale* (SI) and is also recommended by the Italian Consensus document [33,105]. Some laboratories (8% in the CAP survey) did not even know the type of measurement unit they were using [9,99,116].

The use of different units (FEU or DDU, as well as the different measure units) is challenging for clinicians, causes confusing and potentially leads to the misclassification or misdiagnosis of patients [3,4,105]. The fact that 33% of 1500 US laboratories changed their D-dimer units from those recommended by the manufacturer made the situation even more complicated [99]. These variations add additional confusion to diagnostic reasoning, especially when neighbor laboratories use other units [99]. Clinical laboratories should communicate any change in D-dimer measurement units to their stakeholders in a timely manner. Even more importantly, reaching a universal agreement on a uniform measure unit is a compelling need to achieve harmonization of D-dimer immunoassays [99,103].

Age-specific cutoffs

D-dimer values typically increase in parallel with aging, thus leading to a high proportion of elderly patients with D-dimer levels higher than the conventional cutoff set at $500 \mu\text{g/L}$ FEU [1,4,33,99,103,105,117,118]. The use of age-adjusted cutoffs (i.e. [age-adjusted cutoff, $\mu\text{g/L}$ FEU] = [age in years] $\times$ 10) is now universally recommended for reporting results of D-dimer testing (see below) [33,118–120]. Use of these cutoffs would enable a substantial increase in the positive predictive value

(PPV) without significantly impairing the negative predictive value (NPV), and ultimately improve the clinical usefulness of D-dimer measurement in elderly patients (i.e. aged 50 years or older) with low clinical probability [12,33,118–120]. The test most validated for this purpose is the Vidas[®] assay (bioMérieux) [120]. However, the age-adjusted strategy seems to be less useful with the STA-Liatest[®] assay [121]. The analysis of different populations with different immunoassays may explain, at least in part, such discordances. This emphasizes that the results obtained with specific D-dimer assays cannot be extrapolated from one test to another [121].

However, an international survey on the reporting of D-dimer test results showed that published recommendations on the use of age-adjusted cutoffs have not been broadly implemented (i.e. less than 10% of clinical laboratories have implemented this approach) [105].

Along with the 14 combinations of D-dimer units, the use of age-adjusted cutoffs further complicates clinical decision making – there are nearly 30 different possibilities for reporting D-dimer test results [105]. Worldwide efforts are needed to standardize D-dimer test result reporting, to pave the way for widespread implementation of age-adjusted cutoffs [106]. Specific cutoffs to be used for pregnant women and children have also been proposed [106,122,123]. However, as D-dimer increases physiologically throughout pregnancy, the identification of an accurate D-dimer diagnostic threshold will be challenging [124]. The Diagnosis of PE in Pregnancy (DiPEP) biomarker study recently showed that no biomarker has diagnostic efficiency for diagnosing VTE in pregnancy or in the puerperium [125].

Turnaround time

The TAT is a crucial aspect in D-dimer reporting because this test is mostly used in urgent clinical conditions. The Italian consensus document has recently set an overall TAT <1 h, which seems suitable for managing the vast majority of urgent test requests [33].

The use of D-dimer immunoassays characterized by a wide linearity range (i.e. up to $5000 \mu\text{g/L}$) without performing additional external dilutions, along with the use of faster and validated centrifugation processes, PTS, or reliable POC analyzers present valuable approaches to reduce the TAT [3,33,39,43–45,92,107]. Manual ELISAs do not meet the recommended TAT, whilst recent central laboratory D-dimer assays, with analytical process times ranging between 15–40 min, are more likely to fulfill the 1-h criterion (Table 3).

Notably, owing to a half-life of 6–8 h for D-dimer, repeating the test within this time interval has no clinical

basis [33]. In a European study, 81% of participants stated that they measure D-dimer 24 h per day [98].

Clinical applications

More than 30 years ago, D-dimer was found to be a reliable biomarker of both coagulation activation and fibrin digestion [8]. Since then, D-dimer has been used in several clinical conditions. Most widely, D-dimer is used for ruling out VTE in patients with low-medium clinical probability, where D-dimer is now considered the biochemical gold standard [1,4,33,98,126]. D-dimer testing is also useful for assessing the risk of recurrent thrombosis and for guiding anticoagulant therapy, for diagnosing and monitoring DIC, for excluding acute aortic dissection (AAD), and for predicting and managing thrombotic complications in patients with severe infections and sepsis [1,4,127]. Other applications have been studied (e.g. prognostication of peripheral artery disease, identification of vaso-occlusive crisis in sickle cell disease, screening of intracardiac thrombus, prediction of VTE in sleep apnea), but no clinical validation has so far been achieved for these conditions [4,128,129].

Below, we discuss the role of D-dimer for exclusion of VTE, for prediction of VTE recurrence, for prediction of thrombosis in medically hospitalized patients, and for diagnosis and monitoring of DIC.

We briefly present here some indications discussed in the literature, namely cerebral venous thrombosis (CVT), AAD, and acute mesenteric ischemia (AMI).

Some indications for measuring D-dimer

Cerebral venous thrombosis

According to the statement of the American Heart Association (AHA) and of the American Stroke Association (ASA) guidelines, a normal D-dimer level obtained with a sensitive assay may help to identify patients with a low probability of CVT. However, the class and level of evidence are low (Class IIb; Level of Evidence B) and the use of D-dimer is considered superfluous when the clinical suspicion of CVT is high [130]. Because of the variability of presentation, the diagnosis of CVT is rarely confirmed by neuroimaging investigations, which remains, however, the gold standard for the diagnosis of CVT [131]. Pretest clinical probability scores validated to assist clinician in the diagnosis of CVT are also lacking [131]. The meta-analysis of Dentali et al. (14 studies included; 1134 patients) found a mean sensitivity and specificity for D-dimer of 93.9% and 89.7%, respectively, in patients with suspected CVT. The risk of false-negative D-dimer results included longer duration of symptoms, limited sinus involvement and isolated headache [131].

Future prospective studies are needed to confirm these findings. The recent meta-analysis performed by Alons et al. (eight studies included; 636 patients) showed that D-dimer has a high NPV in low-risk patients with isolated headache for excluding CVT [132]. Low-risk patients were defined by normal neurological examination, normal standard head CT (computed tomography), and absence of risk factors such as puerperium and pregnancy. Sensitivity, specificity, PPV, and NPV for diagnosing CVT were 97.8%, 84.9%, 33.1%, and 99.8%, respectively [132]. Normal D-dimer levels may, therefore, reduce unnecessary neuroimaging investigations.

Acute aortic dissection

Symptoms of AAD typically include back and/or abdominal pain, acute onset of tearing chest, asymmetric blood pressure and widened mediastinum on chest x-ray [133]. However, in many patients, symptoms are nonspecific and a missed diagnosis may be fatal [133]. The diagnosis of AAD currently relies on imaging techniques (i.e. magnetic resonance imaging, echocardiography, contrast-enhanced CT) (Class I; Level of evidence B) [133,134]. The potential usefulness of D-dimer lies in the possibility to rule out AAD in patients with low clinical probability because D-dimer has been reported to be persistently elevated in AAD [4,133,135]. D-dimer testing is also rapid, economical, and more accessible than imaging investigations [136]. The meta-analysis of Cui et al. (five studies included; 743 subjects) found a sensitivity and specificity for the diagnostic of AAD of 94.5% and 69.1%, respectively [133]. Similar performances were obtained in the meta-analysis of Watanabe et al. (12 studies included; 2827 subjects; sensitivity and specificity of 95.2% and 60.4%, respectively) [137]. D-dimer levels can, therefore, be used to rule out AAD in patients with low likelihood of the disease [133,137]. According to the American College of Cardiology Foundation (ACCF) and to the AHA guideline, D-dimer cannot be used to rule out the disease in high-risk patients and they do not recommend D-dimer screening for all patients being evaluated for aortic dissection (AD) [134]. The sensitivity and the failure rate of D-dimer in patients with AD detection risk score (ADD-RS) 0, ≤ 1 , and 1 (high risk of AD) were 100% and 0%, 98.7% and 0.8%, and 97.5% and 4.2%, respectively [138]. The ADD-RS allows standardized assessment of pretest probability for acute aortic syndromes [139]. The combined ADD-RS (0 or ≤ 1) with negative D-dimer was superior to D-dimer alone in AAD diagnosis [138,140,141]. However, the accuracy of D-dimer was found to be lower (4% failure rate) in patients at high risk (ADD-RS of 1); therefore D-dimer was not

acceptable to rule out a diagnosis of AD in such settings [141]. In conclusion, the absence of ADD risk markers combined with a negative D-dimer result argues strongly against the diagnosis of AD [138,141]. Integration of ADD-RS with D-dimer may be considered to standardize the diagnosis of AD [141].

Acute mesenteric ischemia

The interest in D-dimer in the diagnosis of AMI has also been acknowledged in the literature [142–145]. Four etiologies of AMI exist: arterial thrombosis (15–20%), venous thrombosis (5%), arterial embolism (50%), and non-occlusive mesenteric ischemia (20–30%) [146,147]. Failure to recognize AMI is responsible for a high mortality rate [147]. The sensitivity and specificity of spiral CT for the detection of AMI have been determined to be 93.3% and 95.9%, respectively [148]. According to the guideline of the European Society for Trauma and Emergency Surgery (ESTES), published in 2016, D-dimer does not discriminate patients with AMI from non-AMI patients [147]. Moreover, no correlation between serum D-dimer levels and the severity of AMI was observed [143]. However, the guideline cited only the paper of Chiu et al. (prospective study on 67 patients), which was published in 2009 [143,147]. In 2017, a meta-analysis performed by Sun et al. included more studies (12 studies; 1300 patients), eight of which concerned AMI. The four other studies dealt with acute strangulated intestinal obstruction, acute intestinal necrosis and mixed types of acute intestinal ischemia, which could be gathered into the term “acute intestinal ischemia”; the calculated sensitivity and specificity were 94% and 50%, respectively [149]. D-dimer measurement could hence be useful for identifying patients with AMI. Further validations of these results in large multi-center clinical studies are needed.

Exclusion of the diagnosis of venous thromboembolism

Although D-dimer levels are increased in almost all cases of acute VTE, any other condition that increases

fibrin production or breakdown may also generate increased D-dimer values [2]. These conditions, listed in Table 4, include infections, cancer, chronic inflammation, aging, pregnancy, recent surgery, and trauma. A large retrospective study of 1647 patients showed that the most common causes of positive D-dimer results were related to infections, followed by VTE, syncope, heart failure, trauma, and cancer [126]. This is of particular concern for hospitalized patients, because D-dimer may be increased for reasons other than VTE and DIC [1]. Brotman et al. found that only 22% of hospitalized patients had normal D-dimer levels [150]. Therefore, D-dimer testing should not be considered specific for VTE [1,2,4,33], and it should be used mainly to exclude VTE because low D-dimer values reflect the lack of significant (ongoing) activation of blood coagulation [1,151].

Epidemiology of VTE

VTE, including both DVT and PE, is a life-threatening condition, with an incidence ranging between 104 to 183 per 100,000 person-years in Europe [152]. The incidence rate is generally higher in men after 45 years, as well as in women in the childbearing years [152]. The prevalence of VTE has been estimated at 422 cases per 100,000 individuals with isolated DVT (69.9%), PE (23.7%) or both (6.4%) [117]. The frequency of VTE cases appears to be slightly higher in women than in men in the United States [117]. The incidence of VTE is also known to increase sharply with age [153].

Of note, during the last decades, the overall incidence of VTE decreased significantly [154]. The wider introduction of thromboprophylaxis in hospitalized patients may, at least in part, explain this result. Specifically, the incidence of isolated DVT also decreased while the incidence of isolated PE increased [154]. This discordance may be explained by the lack of major improvement in the accuracy of ultrasonographic diagnosis of DVT whereas improved resolution imaging has been used for the diagnosis of PE (i.e. multi-detector CT scan) [154].

Table 4. Causes of D-dimer elevation [1,126,127,164].

Acute respiratory distress syndrome	Disseminated intravascular coagulation	Pancreatitis
Advance age	Heart failure	Post transplantation complications
Alzheimer	HELLP syndrome	Pregnancy or puerperium
Aneurysm	Hemolysis (falciform anemia)	Recent surgery
Aortic dissection	Hemorrhage	Renal disease
Arthritis	Hospitalization	Severe chronic urticaria
Atrial fibrillation	Inflammatory bowel disease	Thrombolytic therapy
Burns	Ischemic cardiopathy	Trauma
Cancer	Liver disease	Venous or arterial thrombosis
Chronic inflammation	Localized or systemic Infection	
Disability	Neonatal period	

Several factors have been associated with an increased risk of VTE, including major surgery, hospitalization for acute medical illness (e.g. heart failure, diabetes or pneumonia), trauma/fracture, thrombophilia, obesity, active cancer, pregnancy or postpartum, extended periods of immobility, and oral contraceptives [152,153,155]. Overall, the survival after PE is predictably worse than after DVT alone. The survival rate for DVT at 3 months and 5 years is 91.9% and 72.7%, respectively, whilst for PE, it is 62.8% and 47.5%, respectively [152,156].

VTE is a major healthcare issue worldwide, with an estimated incidence of 3.0–3.3 cases per 100 hospitalizations per year and a rate of hospital readmission of 5–14% [152,157–159]. Compared to hospitalized controls, adjusted mean predicted costs were 1.5–2.5 times higher for patient with VTE. Cost differences were greatest within the first 3 months [160,161]. The total annual health care cost for VTE in 2007 ranged from \$7594 to \$16,644 per patient in the US [158]. The VTE-attributable cost has been estimated to be up to \$10 billion [162]. It has also been estimated that the number of cases of VTE among adults will at least double by 2050 (0.95 million in 2006 to 1.82 million predicted in 2050) [117].

D-dimer assay performance

According to the US Food and Drug Administration (FDA), the specifications for sensitivity and NPV of a D-dimer assay applied to rule out VTE should be $\geq 95\%$ (with lower limit of CI $\geq 90\%$) and $\geq 97\%$ (with lower limit of CI $\geq 95\%$), respectively [163]. The CLSI provides slightly more stringent recommendations for sensitivity and NPV, entailing $\geq 97\%$ (with lower limit of CI $\geq 90\%$) and $\geq 98\%$ (with lower limit of CI $\geq 95\%$), respectively [163]. Based on these requirements, D-dimer assays can be classified into two main categories: (1) those with high sensitivity ($>95\%$) and low specificity ($<40\%$), and (2) those with moderate sensitivity (80–94%) and high specificity (up to 70%) [164]. Quantitative laboratory assays are more likely to have a higher NPV compared to semi-quantitative and qualitative assays [3]. The meta-analysis published by Di Nisio et al., which included 113 studies, showed that ELFAs (e.g. Vidas[®], Stratus DS[®] (Siemens, USA), AxSYM[®]), microplate ELISAs (e.g. Asserachrom[®], Enzygnost[®]), and latex quantitative immunoassays (second generation of latex-based assays [2]) (e.g. Tina-quant[®], STA-Liatest[®]) had higher sensitivity (median, $\geq 95\%$) but lower specificity (median, $\pm 50\%$) [165]. They also found that whole blood agglutination, latex semi-quantitative or qualitative assays had low sensitivity (median, 84%) and high specificity

(median, 76%) [165]. These results were also confirmed in the meta-analysis of Stein et al., which included 78 studies [151]. Similarly, Roy et al. found that ELFA assays exhibit the highest exclusion value (the lowest negative likelihood ratio, see below) and whole blood agglutination assays, the lowest [166]. The specificity of whole blood POC D-dimer assays is overall higher than central laboratory D-dimer assays (71% and 69% for DVT and PE, respectively), meaning that the number of patients in whom VTE could be excluded increases [91]. However, the sensitivity is lower (83% and 87% for DVT and PE, respectively) as compared to central laboratory D-dimer assays [165]. More recent studies also found optimal diagnostic performance for current laboratory D-dimer immunoassays (e.g. STA-Liatest[®], HemosIL HS 500[®], HemosIL AcuStar[®], Innovance[®] (Siemens, USA), Tina-Quant[®]) [38,48,84,112,114,167,168]. Not all POC D-dimer assays are appropriate to be used to exclude VTE [107]. Some POC assays have been approved only as an aid to diagnosis and still require additional tests to rule-out VTE reliably [107]. It is also important to know if the FDA has cleared the assay for exclusion of VTE [107]. Up to now, few POC D-dimer assays meet the performance characteristics recommended by CLSI [95,107]. Quantitative POC assays have been found to be more reliable than qualitative POC assays, especially in patients with low pretest probability of VTE [91,107,169]. Taken together, GPs should weigh the benefit of rapid POC D-dimer assay results over the drawback of overall lower performance [91].

Notably, the NPV is directly influenced by the prevalence of a disease [38]. For example, the prevalence of DVT ranges between 10–50% [38]. An accurate interpretation of D-dimer immunoassay performance is hence needed. The negative likelihood ratio may be the most accurate criterion to assess the diagnostic performance of a D-dimer test [170,171]. The negative likelihood ratio represents the extent of change in the odds or the probability of the disease after a test result [172]. The number needed to test (NNT) index, or the number of patients in whom D-dimer must be measured to rule out one VTE, may also be used to compare the diagnostic yield of commercial D-dimer immunoassays [168,173,174].

Clinical prediction rules (CPR)

It has been demonstrated extensively that assessment by means of a clinical score should be combined with D-dimer testing to improve the diagnostic specificity of the D-dimer result and enhance the diagnostic accuracy [2,4,33,120,153,175]. The aim of clinical probability assessment is to deal with the prevalence dependence

Table 5. Clinical prediction rules for PE (adapted from Righini et al. [174]).

Wells score		Geneva score		Revised Geneva score	
Items	Score	Items	Score	Items	Score
Previous PE or DVT	1.5	Previous PE or DVT	2	Age >65	1
Heart Rate >100	1.5	Heart Rate >100	1	Previous PE or DVT	3
Recent surgery or immobilization	1.5	Recent Surgery	3	Surgery or fracture within 1 month	2
Clinical signs of DVT	3	Age		Active malignancy	2
Alternative diagnosis less likely than PE	3	60–79	1	Unilateral lower limb pain	3
Hemoptysis	1	≥80	2	Hemoptysis	2
Cancer	1	Arterial blood gases		Heart rate	
		CO ₂ (kPa)		75–94	3
		<4.8	2	≥95	5
		4.8–5.19	1	Pain on lower limb deep vein palpation and unilateral edema	4
		O ₂ (kPa)			
		<6.5	4		
		6.5–7.99	3		
		8–9.49	2		
		9.5–10.99	1		
		Chest X-ray			
		Atelectasis	1		
		Elevated hemidiaphragm	1		
Clinical probability	Clinical probability	Clinical probability			
Low	<2	Low	0–4	Low	0–3
Intermediate	2–6	Intermediate	5–8	Intermediate	4–10
High	>6	High	≥9	High	≥11
Dichotomized					
PE unlikely	≤4				
PE likely	>4				

for accurate interpretation of diagnostic tests of D-dimer assays by identifying subgroups of patients according to different levels of VTE prevalence: low (<10%), moderate (around 30%), high (>50%), or "unlikely" ($\pm 10\%$) and "likely" ($\pm 35\%$) [155]. There are additional aspects that support the use of D-dimer in combination with CPR, such as the fact that some patients with thrombosis may present a false-negative D-dimer concentration due to a hypofibrinolytic state, small thrombi (i.e. distal DVT or isolated subsegmental PE), anticoagulant therapy or D-dimer testing performed too early or too late after the thrombosis [1,2,4,33,48].

CPR have been derived from clinical studies and include a number of signs, symptoms, and risk factors of VTE [2]. The use of two-level scores is typically preferred, because it is straightforward and because three-level scores depend on the sensitivity of the D-dimer assay (which is rarely known by the clinicians) [153,175]. The Well's and Geneva scores (modified or not) are the two CPRs mainly used for PE, whilst the Well's score is especially used for DVT [155,169,174,176]. These scores for PE are shown in Table 5. The diagnostic performances of these CPR were found to be almost comparable [119,177–179]. Other CPR have also been validated and may be used routinely (e.g. Empiric, Charlotte or Simplified Geneva) [1,177,180]. A good CPR should not contain subjective variables, and it should be accurate, reproducible, easy to remember, and offer a standardized approach compared to clinical assessment (especially for

inexperienced physicians) [1,155]. However, pretest probability calculation may still be subjective [181,182]. These algorithms are considered clinically validated when the number of thromboembolic event after 3 months is as high as 1–2% compared to a gold standard test (pulmonary angiography for PE and venography for DVT) [183].

Different algorithms combining CPR and D-dimer measurement have been proposed and endorsed by international guidelines [33,155,164,184]. A "low" or "unlikely" probability is generally followed by D-dimer testing. VTE can be ruled out when D-dimer values are below the method-specific cutoff. In these conditions, the use of time-consuming, expensive and less safe (i.e. due to radiation exposure and/or injection of contrast products) tests would be unnecessary [2]. However, when the clinical probability is "high or likely" or "moderate/high", D-dimer testing is not useful and radiologic imaging should be used [2,185]. A positive D-dimer test result may also trigger imaging tests [2]. For DVT, this would entail ultrasonography [183] or Doppler flow studies [183]. For PE, the most used techniques include CT pulmonary angiogram (CTPA) [183] or contrast-enhanced or unenhanced magnetic resonance imaging, especially when CTPA is unadvisable [183]. Some studies based on specific algorithms concluded that a negative D-dimer result in patients with a "moderate" score could also safely rule out a VTE episode, thus enhancing the usefulness of D-dimer testing and decreasing the use of imaging testing [155,180]. Several other studies, including a meta-analysis of 1660

patients, showed that a normal D-dimer concentration along with a low CPR score had a NPV of up to 95% [3,116,119,186].

However, in current practice, such clinical algorithms are not always followed [181,182]. In a recent survey that included 487 clinicians from six countries, only 70.3% used pretest probability scores, and 10% excluded or confirmed DVT based only on D-dimer test results. Moreover, a significant number of clinicians still order D-dimer testing in patients with a high probability of VTE, whilst others order imaging testing in cases of low pretest probability. This waste of resources may lead to a negative impact on health care resources [181,182]. Moreover, in a large study that included 117 centers and 1529 patients with suspected PE, the diagnostic management was inappropriate in 662 patients (43%): in 36 of 429 (8%) patients with confirmed PE and in 626 of 110 (57%) patients in whom PE was ruled out [187]. One of the independent risk factors for inappropriate management was the lack of written algorithms that included a clinical probability score. Major efforts are hence needed to harmonize the use of the different CPR algorithms.

A recent prospective study based on 808 patients found that a low D-dimer level ($<750 \mu\text{g/L}$ FEU) measured with the MDA D-dimer assay (quantitative latex agglutination assay, bioMérieux) could be used alone for ruling out PE, without the need of a CPR score [188]. The theoretical advantage of limiting the use of CPR is a decrease in the cost attributable to imaging techniques (e.g. CTPA and ventilation-perfusion lung scanning) and prevention of radiation exposure [188]. However, the MDA assay is no longer being marketed. Moreover, the number of patients with a high pretest probability was quite low in that study [188]. On the other hand, D-dimer testing has limited clinical value in patients with a high CPR score in whom the prevalence of VTE is high (i.e. low NPV) [164]. Additional studies with other D-dimer assays are needed to confirm these findings. Nevertheless, the use of a D-dimer assay as a stand-alone test has not been widely endorsed by most guidelines and consensus documents [1,2,4,33,38,127,164,175,183,189]. This is understandable because the high percentage of deaths in the untreated patients is the major risk of underdiagnosing VTE [183,185].

Age-adjusted cutoff

As earlier discussed, D-dimer values tend to increase with age, thus leading to a greater proportion of older patients ($\approx 60\%$) with D-dimer values higher than the conventional cutoff of $500 \mu\text{g/L}$ FEU [1,4,33,105,117,118,190]. The specificity of D-dimer for VTE in elderly patients (0–18%,

≥ 80 years) is lower compared to younger patients (49–67%, <50 years). Elderly patients are therefore less likely to have negative D-dimer results when VTE is absent [12,118], and a high number of these patients with a low clinical score may undergo unnecessary imaging testing because of high D-dimer values [118]. This is because consolidated evidence shows that VTE prevalence increases with age [12,119,152,174].

Age-adjusted cutoffs have subsequently been validated and recommended because their application increases diagnostic specificity without significantly impairing the NPV [33,119,173,191–193]. The calculation of the age-adjusted cutoff is based on the following equation: "age-adjusted cutoff, $\mu\text{g/L}$ FEU = age in years $\times 10$ [119]". Thus, the age-adjusted cutoff of a 70-year-old patient is $700 \mu\text{g/L}$ FEU (rather than $500 \mu\text{g/L}$ FEU). The multiplication factor of "10" has been slightly rounded to facilitate practicability and clinical usefulness [119].

The ADJUST-PE study, which included over 3000 patients, prospectively demonstrated that the combination of age-adjusted cutoff and a clinical score (Wells or Revised Geneva score) was associated with a higher proportion of patients in whom PE could be safely ruled out without additional CTPA [120,194]. Many other studies have demonstrated that the number of patients in whom VTE could be ruled out increased remarkably after applying age-adjusted cutoffs, and this approach has also been found to be cost-effective because a lower number of patients with low clinical score will need to undergo expensive imaging tests [120,168,174,183,194].

A recent meta-analysis evaluated the clinical effectiveness of age-adjusted cutoffs in 12497 older patients and concluded that the sensitivity was $>97\%$ in all age categories [118]. The gain in specificity (up to 23% better specificity compared to a fixed cutoff) was especially favorable for patients aged 75 years or older [118,120]. This strategy has been validated for both PE and DVT with many commercial D-dimer immunoassays [38,119,120,168,191]. Specifically, Oude Elferink et al. found that the NPV of 9 D-dimer immunoassays was always $>97\%$ [38]. Another study found a NPV $>99\%$ in using five different D-dimer immunoassays [168].

More recently, Farm et al. reported that 4 D-dimer assays display sensitivity $>95\%$ and NPV $>97\%$, whilst only one method (MediRox[®], Sweden) displayed slightly lower sensitivity (94%) [191]. Although age-adjusted cutoffs seem straightforward in most cases, caution should be used with non- or insufficiently validated D-dimer immunoassays [106]. More recently, the use of age-adjusted cutoffs for ruling out PE has also been endorsed by the European Society of Cardiology (ESC) [184].

Interestingly, a recent study demonstrated that age-adjusted cutoffs may also be useful in younger patients. Prochaska et al. [192] performed a clinical study that included 500 patients with suspected DVT, and reported that a D-dimer cutoff of 500 µg/L FEU had lower sensitivity in patients <60 years (76.1%) compared to older patients (92.9%). A lower sensitivity was found especially in women, and in patients with unprovoked DVT, with low thrombotic burden, and with distal DVT [192]. However, the use of a lower cutoff in younger patients (250 µg/L FEU) was effective in increasing the sensitivity of this test remarkably, up to 92.4%. Nevertheless, further studies will be needed to confirm the results of this preliminary observational study.

Clinical probability-adjusted cutoffs

The use of clinical probability-adjusted cutoffs has also been proposed [2,195–197]. Linkins et al. used a low (200 µg/L FEU), an intermediate or “conventional” (500 µg/L FEU) and a high (2000 µg/L FEU) D-dimer cutoff (MDA D-dimer assay) for high, moderate, and low pretest probability patients, respectively [196]. Overall, 571 patients were studied retrospectively. Although the authors failed to find different sensitivities and NPVs across the groups for diagnosing VTE, the use of clinical probability-adjusted cutoffs increased the D-dimer test specificity compared to the conventional approach (i.e. 44.7% versus 60.4%) [196]. Another study published by Kabrhel et al. in 7940 patients with suspected PE and based on six different D-dimer assays reported similar data [195]. The use of higher (conventional $\times 2$), conventional (identical), and lower (conventional/2) cutoffs for low, intermediate, and high pretest probability patients displayed a higher specificity (75% versus 58%) and a similar NPV (99.1% versus 99.5%) [195]. Kline et al. also suggested using a higher D-dimer cutoff (1000 µg/L instead of 500 µg/L) in younger patients (i.e. aged <50 years) with a low clinical score for PE (clinical probability-adjusted strategy). The authors found a significant decrease in the need for performing unnecessary pulmonary angiograms but also observed a slight decrease in diagnosis of isolated subsegmental PE [198]. Likewise, Linkins et al. reported that modifying D-dimer cutoffs according to the clinical pretest probability of DVT (<500 µg/L and <1000 µg/L for moderate and low pretest probability, respectively) may be associated with a reduced need to perform ultrasonography, without impairing the NPV of the test [185].

Recently, van der Hulle et al. prospectively validated a diagnostic algorithm incorporating the YEARS clinical decision rule (CDR) (clinical signs of DVT, hemoptysis,

PE as the most likely diagnosis) and various D-dimer cutoffs in 3465 patients with suspected PE. Among the 2946 (85%) patients in whom PE was ruled out (no YEARS CDR and D-dimer <1000 µg/L or at least one YEARS CDR and D-dimer <500 µg/L), the rate of VTE during a 3-month follow-up period was only 0.61% (95% CI: 0.36%–0.96%) [180]. Furthermore, a 14% decrease in CTPA was observed [180]. Based on this latter study, van der Pol et al. found that the age-adjusted cutoff added no value to the YEARS algorithm [199].

Takach Lapner et al. recently compared both the age-adjusted and the clinical probability-adjusted strategies [197]. For that purpose, 1649 patients with low or moderate clinical probability after assessment for a first VTE event were studied retrospectively. The authors found that the NPV was similar in both groups (99.7% versus 99.6%), but also observed a significantly higher rate of negative D-dimer test results (STA-Liatest®) in the clinical probability-adjusted strategy compared to the age-adjusted strategy (56.1% versus 50.9%) [197]. According to these recent findings, it can be concluded that the clinical probability-adjusted cutoffs may be more efficient than the age-adjusted cutoffs for limiting the number of imaging tests. However, this conclusion needs to be supported by further prospective studies using different D-dimer immunoassays.

VTE diagnosis in specific populations

In patients with recent surgery, D-dimers are usually increased and the diagnosis of VTE is hence challenging. No CPR is validated and only the global clinical judgment of the clinician can be used to assess CPR [200]. However, clinical judgment is criticized for its lack of standardization and the difficulty in teaching. Therefore, clinicians perform imaging whenever VTE is suspected. Nevertheless, Penaloza et al. developed a score composed of 9 variables to assess the usefulness of D-dimer measurement, using data from 4537 patients that had high sensitive D-dimer testing (Vidas® for 2079 patients, STA-Liatest® for 1783 patients, and MDA® for 675 patients) [201]. These 9 variables were independently associated with a risk of false positive D-dimer value (sex female: +1, age 65–84 years: +4, age ≥ 85 years: +8, heart rate ≥ 95 /min: +1, pulse oxygen saturation <95%: +2, temperature ≥ 38.5 °C: +3, personal history of VTE: +1, surgery under general anesthesia within 4 weeks: +2, active malignancy: +3, pregnancy or postpartum within 4 weeks: +4). In patients with non-high CPR and a relevance score ≤ 8 , at least 10% had a quantitative D-dimer assay below 500 µg/L, meaning that, in patients with recent surgery, D-dimer measurement may still be useful, at least if

they have no other risk factors for false positive results. Indeed, a negative D-dimer result keeps its sensitivity and NPV [200].

Similarly, D-dimer levels increased physiologically in pregnancy and in the postpartum period. In a study that included 1343 pregnant women with D-dimer measurements using a turbidimetric method (STA-Liatest®), the rate of D-dimer test results below the usual cutoff (500 µg/L) in healthy pregnant women was 85%, 29% and 4.1% during the first, the second and the third trimesters, respectively [202]. Similar results were observed with other assays (MDA® and HemosIL HS assay®) [203,204]. In the post-partum period, D-dimer returns to normal levels around the 6th week. In case of suspicion of PE, because imaging tests may expose the mother and the fetus to radiation, the ability to rule-out PE with non-radiologic tests is crucial [107]. To improve the usefulness of D-dimer results, specific cutoffs have been proposed, but not yet validated, in a large prospective study [124]. However, once again, a negative D-dimer result keeps its sensitivity and NPV [107]. Two large studies are currently ongoing on this topic and may give important information on the best algorithm to exclude VTE in pregnant women [205]. Currently and in order to avoid imaging tests, D-dimer measurement is still recommended in pregnant women with suspicion of VTE and a non-high CPR, at least in the suspicion of PE [206].

In patients known to have cancer, the relevance of D-dimer in the diagnosis of VTE is decreased [107]. The prevalence of VTE is increased (up to 20% of cancer patients develop VTE) and the NPV is therefore reduced [107,207]. A large meta-analysis of 10002 patients showed that the prevalence of both a low Wells score and a negative D-dimer value among patients with cancer was only 9% [208]. It has also been shown that 88–94% of patients with malignancy require additional tests to rule out VTE [209]. However, the D-dimer measurement as included in a score may be useful to identify cancer patients at low or high risk of VTE [210].

The exclusion of VTE in patients with renal disease is also challenging, given that D-dimer fragments are eliminated mainly by renal clearance and reticuloendothelial system catabolism [11]. Accordingly, Robert-Ebadi et al. showed that D-dimer concentrations increased according to renal status, and the proportion of negative D-dimer results decreased from 46% to 11% with normal (≥ 90 mL/min) and moderate renal failure (30–59 mL/min), respectively [211]. In addition, Lindner et al. showed that 100% of patients with low eGFR (< 30 mL/min; $n = 29$) had a positive D-dimer value [212]. The specificity will therefore proportionally

decrease according to the renal failure status [212,213]. The possibility of using renal function-adjusted cutoffs should be verified in dedicated studies [212].

Prediction of recurrence of VTE

The risk of recurrence during a 1-year follow up after a first VTE episode is higher in men (9.5%) than in women (5.3%), and tends also to increase over time (9.1% and 19.7%, for men and women after 3 years, respectively) [214]. Interestingly, a high D-dimer value has been associated with increased risk of recurrent VTE [127]. It has hence been suggested that the duration of anticoagulant therapy in VTE patients may be combined with results of D-dimer testing [1,215]. In the PROLONG study, the D-dimer value was increased 1 month after cessation of vitamin K antagonist (VKA) in 608 VTE anticoagulated patients who were treated for at least 3 months [215]. Patients with normal D-dimer values measured with the Clearview Simplify® D-dimer assay (whole blood qualitative assay) did not resume anticoagulation, whilst those with abnormal D-dimer values were randomly assigned to receive or not receive a new anticoagulant therapy. An abnormal D-dimer concentration was observed in $\approx 37\%$ of patients. Among these, a new VTE episode occurred in 15% of patients who stopped VKA, whilst VTE events occurred only in 2.9% of those receiving VKA therapy (adjusted hazard ratio [HR] 4.26; 95% CI 1.23–14.6; $p = .02$). Legnani et al. reanalyzed the data of the PROLONG study to assess D-dimer cutoff values for several quantitative assays (Vidas®, Innovance®, HemosIL®, and STA-Liatest®) [216]. Based on this post-hoc analysis, they determine age-specific cutoffs and identified a higher threshold in older patients (i.e. aged 70 years or older) [216]. The use of a strategy based on method-specific cutoff values, adjusted for age and gender, is therefore advisable for more accurate prediction of VTE recurrence [217,218]. The use of quantitative methods rather than qualitative tests results (i.e., positive or negative) may be more suitable to achieve this goal [216–218]. The meta-analysis of Douketis et al., which included 1818 patients from seven prospective studies, found that the D-dimer value was a significant predictor of VTE (HR 2.59; 95% CI 1.90–3.52; $p < .001$) [214]. Overall, a D-dimer value higher than the diagnostic cutoff after 3-months of anticoagulant therapy in patients with a first unprovoked VTE event was associated with double the risk of recurrence compared to patients with lower D-dimer values [2,219]. The risk of VTE recurrence was similar among young and elderly patients. More recently, Nagler et al. showed that increased D-dimer values

(Vidas® and Innovance®) observed after discontinuing anticoagulant therapy after 1 month (HR 3.3; 95% CI 1.8–6.1), male sex (HR 2.8; 95% CI 1.5–5.1), and use of oral contraceptives (HR 0.1; 95% CI 0.0–0.9), were significant predictors of recurrent DVT in 479 patients diagnosed with proximal DVT [220]. These authors also found that high factor VIII values (HR 2.2; 95% CI 1.2–4.0) could be predictors of recurrent DVT [220]. The distinction between provoked and unprovoked VTE is also important. Patients with provoked VTE (after surgery, lower limb trauma, and orthopedic immobilization or hospitalization for an acute medical illness) have a very low risk of recurrence as compared to patients with unprovoked VTE who have a higher risk of recurrence at two years ($\pm 20\%$) [221].

Several models have been developed to predict the risk of recurrence of unprovoked VTE. In the real-life study of Tosetto et al., performed with different D-dimer assays in nine centers, an increased D-dimer value after stopping anticoagulation (cutoff 500 $\mu\text{g/L}$ with a quantitative assay or “positive” with a qualitative assay), a younger age (<50 years), male sex, and VTE not associated with hormonal therapy in women were identified as the best predictors of VTE recurrence (area under the curve, 0.71) [222]. The diagnostic performance of a model including all these variables had a higher diagnostic efficiency compared to D-dimer testing alone (area under the curve 0.61; $p < .001$). This score may hence be useful to assess whether the anticoagulant treatment should be discontinued or resumed after the usual 3-month period. The DASH score has been retrospectively validated in a recent independent cohort study, but the score appeared to perform better in younger subjects (<65 years old) [223]. To our knowledge, the DASH score has not been validated in an interventional prospective study. Additional algorithms, such as the Vienna and HERDOO2 prediction models, have also been developed [224]. The Vienna prediction model included D-dimer (Asserachrom assay®), sex, and site of index event [224,225]. It has been validated in external populations but, again, not in an interventional study [226]. Conversely, the HERDOO2 model has been recently validated in an interventional study that included 3155 patients from several countries [227]. This model included D-dimer (measured with Vidas® assay; cutoff, 250 $\mu\text{g/L}$), age, body mass index (BMI), and post-thrombotic signs [224,228]. In line with these findings, the Scientific and Standardization Committee of the International Society on Thrombosis and Hemostasis recommends that D-dimer testing should be performed in all patients with clinical suspicion of recurrent VTE

[229]. However, caution has to be used when translating these findings to other D-dimer immunoassays. Moreover, controversy remains about the most appropriate timing of D-dimer monitoring during or after discontinuing anticoagulation. Finally, the accurate definition of “positive” or “negative” D-dimer (i.e. precise cutoff) is still an unmet requirement [2,214].

Prediction of risk of VTE in hospitalized patients

Recently, a baseline D-dimer level measured in 7441 hospitalized patients was found to be independently associated with symptomatic VTE during a period of ~ 3 months (adjusted HR 2.22; 95% CI 1.38–1.58; $p < .001$) [230,231]. Patients hospitalized for acute illness (ischemic stroke, heart failure, respiratory failure, rheumatic disorders, infection) have an increased risk of developing VTE [231–234]. A risk assessment for VTE is hence needed before initiating thromboprophylaxis [231]. The International Medical Prevention Registry on Venous Thromboembolism (IMPROVE) assessment tool has been developed with the aim of risk stratification in hospitalized, medically ill patients. The score encompasses many clinical variables such as previous VTE, known thrombophilia, current lower-limb paralysis, current cancer, immobilization ≥ 7 days, intensive care or coronary care units, and age > 60 years [231]. Because D-dimer is an independent predictor of symptomatic VTE, Gibson et al. combined D-dimer testing with the IMPROVE assessment tool [231]. A D-dimer level $\geq 2 \times$ the upper reference limit (URL) measured with STA-Liatest platform was assigned two more points in the new scoring system (IMPROVEDD). The population of the APEX study was used to validate this approach. The combination of clinical variables already included in the IMPROVE score with D-dimer test results was effective in substantially improving risk discrimination and patients risk reclassification at both 42 and 77 days.

Diagnosis and monitoring of disseminated intravascular coagulation

DIC is a life-threatening condition characterized by persistent activation of the hemostatic system with intravascular thrombin generation, fibrin formation, and increased fibrinolysis [3]. Patients with DIC may present with bleeding, thrombosis, or both [164]. Early recognition of DIC is paramount to initiate the appropriate treatment, which generally entails eliminating the underlying condition (e.g. sepsis, malignancy, trauma or burns, obstetrical diseases, toxins, drugs, immunological disorders, and other inflammatory diseases) [1]. The

International Society of Thrombosis and Haemostasis (ISTH) [235] and the Japanese Ministry of Health and Welfare [236] have proposed two different scoring systems for diagnosis and management of DIC. Comparison of these scoring systems yielded similar diagnostic effectiveness for DIC [1,237]. The laboratory parameters included in these scoring systems are platelet count, fibrinogen, prothrombin time, and FDP [235,236,238]. Although D-dimer is the most common FDP measured in clinical laboratories, some authors suggested that soluble fibrin may be more specific to detect intravascular clot generation. D-dimer is considered positive when its value is ≥ 2 times the URL [1,4]. As for VTE, there is now consolidated evidence that a "normal" D-dimer value excludes a diagnosis of DIC [127]. However, a "positive" D-dimer value cannot be considered specific for DIC, because D-dimer may be increased in a vast array of clinical conditions (Table 4). Repeated D-dimer analysis (i.e. 3 times per day) is also important to monitor DIC [236,238,239]. The combination of D-dimer and fibrin monomer allowed the identification of patients with septic shock with a worse survival [240].

Conclusions

Today, D-dimer has become one of the most commonly requested coagulation tests. D-dimer, considered a biomarker of activation of coagulation and fibrinolysis, is mainly employed for the exclusion of VTE. A negative D-dimer value along with a low clinical probability can safely exclude VTE in suspected individuals. Currently, it is recommended to employ D-dimer assays with a very high sensitivity and NPV to safely exclude VTE ($\geq 95\%$ and $\geq 97\%$, respectively). However, a positive D-dimer value should always trigger imaging tests to confirm VTE, given that a wide array of diseases and conditions are characterized by an increase in D-dimer (i.e. aging, inflammation, cancer, renal failure). More recently, the use of age-adjusted cutoffs to safely rule out VTE has been proposed to substantially increase the PPV without significantly impairing the NPV, so ultimately improving the clinical usefulness of D-dimer measurement in elderly patients with low clinical probability. Even though these adjusted cutoffs have been endorsed in several guidelines, many laboratories are still not using them. Major efforts are therefore needed for implementation of these recommendations. VTE exclusion in specific populations (i.e. pregnancy, heart failure, cancer, perioperative settings) is also challenging. Utilization of POC D-dimer assays as a rapid screening tool is especially attractive for GPs because central laboratory quantitative D-dimer assays are not always

readily available. However, not all POC D-dimer assays are appropriate to be used to exclude VTE, and only POC assays that have been validated in clinical trials and cleared by the FDA should be used.

D-dimer measurement is also useful to predict the risk of VTE recurrence and the risk of VTE in hospitalized patients and to diagnose and monitor DIC. Other emerging indications (i.e. prognostication of peripheral artery disease, screening of intracardiac thrombus) need further investigations.

Pre-analytical requirements are paramount in the hemostasis laboratory. Nevertheless, the literature appears to be quite reassuring regarding several pre-analytical variables that affect D-dimer results (needle bore size, butterfly devices, use of PTS, *in vitro* hemolysis, stability). In some specific settings (i.e. emergency, older patients, delay before centrifugation), D-dimer results may be released to the clinician without impairing patient management.

A major drawback of D-dimer assays is the high variability observed between immunoassays. This variability is explained by the fact that D-dimers comprise a broad mixture of degradation products of cross-linked fibrin, by the utilization of different monoclonal antibodies, by the lack of international certified internal controls or calibrators, and by the use of different units and clinical cutoffs. Further studies are needed to achieve harmonization of D-dimer measurements. Ongoing discussion among manufacturers, scientists and clinicians is essential for achieving this goal.

Acknowledgments

The authors would like to thank Mr. Nicolas Bailly for the editing of the figures and tables.

Disclosure statement

No potential conflict of interest was reported by the authors.

ORCID

Giuseppe Lippi  <http://orcid.org/0000-0001-9523-9054>
François Mullier  <http://orcid.org/0000-0001-6947-6099>

References

- [1] Adam SS, Key NS, Greenberg CS. D-dimer antigen: current concepts and future prospects. *Blood*. 2009; 113:2878–2887.
- [2] Linkins LA, Takach Lapner S. Review of D-dimer testing: good, bad, and ugly. *Int J Lab Hematol*. 2017;39 (Suppl 1):98–103.

- [3] Riley RS, Gilbert AR, Dalton JB, et al. Widely used types and clinical applications of D-dimer assay. *Lab Med.* 2016;47:90–102.
- [4] Thachil J, Lippi G, Favaloro EJ. D-dimer testing: laboratory aspects and current issues. *Methods Mol Biol.* 2017;1646:91–104.
- [5] Blomback B. Fibrinogen and fibrin—proteins with complex roles in hemostasis and thrombosis. *Thromb Res.* 1996;83:1–75.
- [6] Sidelmann JJ, Gram J, Jespersen J, et al. Fibrin clot formation and lysis: basic mechanisms. *Semin Thromb Hemost.* 2000;26:605–618.
- [7] Medved L, Nieuwenhuizen W. Molecular mechanisms of initiation of fibrinolysis by fibrin. *Thromb Haemost.* 2003;89:409–419.
- [8] Gaffney PJ, Lane DA, Kakkar VV, et al. Characterisation of a soluble D dimer-E complex in crosslinked fibrin digests. *Thromb Res.* 1975;7:89–99.
- [9] Reber G, De Moerloose P. Standardization of D-dimer testing. In: Kitchen S, Olson JD, Preston FE, editors. *Quality in laboratory hemostasis and thrombosis*. 2nd ed. Oxford, UK: John Wiley & Sons; 2013. p. 136–146.
- [10] Gaffney PJ, Edgell T, Creighton-Kempsford LJ, et al. Fibrin degradation product (FnDP) assays: analysis of standardization issues and target antigens in plasma. *Br J Haematol.* 1995;90:187–194.
- [11] Lip GY, Lowe GD. Fibrin D-dimer: a useful clinical marker of thrombogenesis? *Clin Sci.* 1995;89:205–214.
- [12] Haase C, Joergensen M, Ellervik C, et al. Age- and sex-dependent reference intervals for D-dimer: evidence for a marked increase by age. *Thromb Res.* 2013;132:676–680.
- [13] Dempfle CE. Validation, calibration, and specificity of quantitative D-dimer assays. *Semin Vasc Med.* 2005;5:315–320.
- [14] ISO. Medical laboratories—Requirements for quality and competence. ISO 15189:2012. Geneva, Switzerland: International Organization for Standardization; 2012.
- [15] Lippi G, Franchini M, Montagnana M, et al. Quality and reliability of routine coagulation testing: can we trust that sample? *Blood Coagul Fibrinolysis.* 2006;17: 513–519.
- [16] Plebani M. Errors in clinical laboratories or errors in laboratory medicine? *Clin Chem Lab Med.* 2006;44: 750–759.
- [17] Plebani M, Carraro P. Mistakes in a stat laboratory: types and frequency. *Clin Chem.* 1997;43:1348–1351.
- [18] Plebani M, Favaloro EJ, Lippi G. Patient safety and quality in laboratory and hemostasis testing: a renewed loop? *Semin Thromb Hemost.* 2012;38: 553–558.
- [19] Preston FE, Lippi G, Favaloro EJ, et al. Quality issues in laboratory haemostasis. *Haemophilia.* 2010;16: 93–99.
- [20] Lippi G, Favaloro E. Causes of errors in medical laboratories. In: Kitchen S, Olson JD, Preston FE, editors. *Quality in laboratory hemostasis and thrombosis*. 2nd ed. Oxford, UK: John Wiley & Sons; 2013. p. 22–31.
- [21] Magnette A, Chatelain M, Chatelain B, et al. Pre-analytical issues in the haemostasis laboratory: guidance for the clinical laboratories. *Thromb J.* 2016;14:49.
- [22] Salvagno GL, Lippi G, Bassi A, et al. Prevalence and type of pre-analytical problems for inpatients samples in coagulation laboratory. *J Eval Clin Pract.* 2008;14:351–353.
- [23] Grecu DS, Vlad DC, Dumitrascu V. Quality indicators in the preanalytical phase of testing in a stat laboratory. *Lab Med.* 2014;45:74–81.
- [24] Dikmen ZG, Pinar A, Akbiyik F. Specimen rejection in laboratory medicine: necessary for patient safety? *Biochem Med.* 2015;25:377–385.
- [25] Adcock D. Sample integrity and preanalytical variables. In: Kitchen S, Olson JD, Preston FE, editors. *Quality in laboratory hemostasis and thrombosis*. 2nd ed. Oxford, UK: John Wiley & Sons; 2013. p. 45–56.
- [26] Fraser CG. *Biological variation: from principles to practice*. Washington, DC: AACC Press; 2001. p. 67–90.
- [27] de Maat MP, van Schie M, Kluft C, et al. Biological variation of hemostasis variables in thrombosis and bleeding: consequences for performance specifications. *Clin Chem.* 2016;62:1639–1646.
- [28] Favaloro EJ, Funk DM, Lippi G. Pre-analytical variables in coagulation testing associated with diagnostic errors in hemostasis. *Lab Med.* 2012;43:1–10.
- [29] Hunt BJ, Parratt R, Cable M, et al. Activation of coagulation and platelets is affected by the hydrophobicity of artificial surfaces. *Blood Coagul Fibrinolysis.* 1997;8:223–231.
- [30] Loeffen R, Kleinegris MC, Loubele ST, et al. Preanalytic variables of thrombin generation: towards a standard procedure and validation of the method. *J Thromb Haemost.* 2012;10:2544–2554.
- [31] Lippi G, Salvagno GL, Montagnana M, et al. Influence of the needle bore size used for collecting venous blood samples on routine clinical chemistry testing. *Clin Chem Lab Med.* 2006;44:1009–1014.
- [32] Lippi G, Salvagno GL, Guidi GC. No influence of a butterfly device on routine coagulation assays and D-dimer measurement. *J Thromb Haemost.* 2005;3: 389–391.
- [33] Lippi G, Cervellin G, Casagrande I, et al. D-dimer testing for suspected venous thromboembolism in the emergency department. Consensus document of AcEMC, CISMEL, SIBioC, and SIMeL. *Clin Chem Lab Med.* 2014;52:621–628.
- [34] CLSI. Collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays. CLSI document H21-A5. 5th ed. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- [35] Leroy-Matheron C, Gouault-Heilmann M. Influence of conditions of blood sampling on coagulation activation markers (prothrombin fragment 1 + 2, thrombin-antithrombin complexes and D-dimers) measurements. *Thromb Res.* 1994;74:399–407.
- [36] Gosselin RC, Janatpour K, Larkin EC, et al. Comparison of samples obtained from 3.2% sodium citrate glass and two 3.2% sodium citrate plastic blood collection tubes used in coagulation testing. *Am J Clin Pathol.* 2004;122:843–848.

- [37] Yavaş S, Ayaz S, Köse SK, et al. Influence of blood collection systems on coagulation tests. *Turk J Haematol.* 2012;29:367–375.
- [38] Oude Elferink RF, Loot AE, Van De Klashorst CG, et al. Clinical evaluation of eight different D-dimer tests for the exclusion of deep venous thrombosis in primary care patients. *Scand J Clin Lab Invest.* 2015; 75:230–238.
- [39] Schutgens RE, Haas FJ, Ruven HJ, et al. No influence of heparin plasma and other (pre)analytic variables on D-dimer determinations. *Clin Chem.* 2002;48: 1611–1613.
- [40] Vukovich TC, Hamwi A, Bieglmayer C. D-dimer testing within the routine clinical chemistry profile. *Clin Chem.* 1998;44:1557–1558.
- [41] Lippi G, Salvagno GL, Rossi L, et al. Analytical performances of the D-dimer assay for the Immulite 2000 automated immunoassay analyser. *Int J Lab Haematol.* 2007;29:415–420.
- [42] Lacroix R, Judicone C, Poncelet P, et al. Impact of pre-analytical parameters on the measurement of circulating microparticles: towards standardization of protocol. *J Thromb Haemost.* 2012;10:437–446.
- [43] Le Quellec S, Paris M, Nougier C, et al. Pre-analytical effects of pneumatic tube system transport on routine haematology and coagulation tests, global coagulation assays and platelet function assays. *Thromb Res.* 2017;153:7–13.
- [44] Wallin O, Soderberg J, Grankvist K, et al. Preanalytical effects of pneumatic tube transport on routine haematology, coagulation parameters, platelet function and global coagulation. *Clin Chem Lab Med.* 2008;46:1443–1449.
- [45] Boudaoud L, Divaret Marie G, Bezeaud PA. Rapid centrifugation for routine coagulation testing. *Ann Biol Clin.* 2006;64:315–317.
- [46] Salvagno GL, Lippi G, Montagnana M, et al. Influence of temperature and time before centrifugation of specimens for routine coagulation testing. *Int J Lab Hematol.* 2009;31:462–467.
- [47] Toulon P, Metge S, Hangard M, et al. Impact of different storage times at room temperature of unspun citrated blood samples on routine coagulation tests results. Results of a bicenter study and review of the literature. *Int J Lab Hem.* 2017;39:458–468.
- [48] Elf JL, Strandberg K, Svensson PJ. Performance of two relatively new quantitative D-dimer assays (Innovance D-dimer and AxsYM D-dimer) for the exclusion of deep vein thrombosis. *Thromb Res.* 2009;124:701–705.
- [49] Lippi G, Avanzini P, Zobbi V, et al. Influence of mechanical hemolysis of blood on two D-dimer immunoassays. *Blood Coagul Fibrinolysis.* 2012;23:461–463.
- [50] Lippi G, Blanckaert N, Bonini P, et al. Haemolysis: an overview of the leading cause of unsuitable specimens in clinical laboratories. *Clin Chem Lab Med.* 2008;46:764–772.
- [51] Lippi G, Plebani M, Favaloro EJ. Interference in coagulation testing: focus on spurious hemolysis, icterus, and lipemia. *Semin Thromb Hemost.* 2012;39: 258–266.
- [52] Lippi G, von Meyer A, Cadamuro J, et al. Blood sample quality. *Diagnosis (Berl).* 2018. DOI: 10.1515/dx-2018-0018.
- [53] Lippi G, Montagnana M, Salvagno GL, et al. Interference of blood cell lysis on routine coagulation testing. *Arch Pathol Lab Med.* 2006;130:181–184.
- [54] D'Angelo G, Villa C, Tamborini A, et al. Evaluation of the main coagulation tests in the presence of hemolysis in healthy subjects and patients on oral anti-coagulant therapy. *Int Jnl Lab Hem.* 2015;37: 819–833.
- [55] Carraro P, Servidio G, Plebani M. Hemolyzed specimens: a reason for rejection or a clinical challenge? *Clin Chem.* 2000;46:306–307.
- [56] Lippi G, Ippolito L, Tondelli MT, et al. Interference from heterophilic antibodies in D-dimer assessment. A case report. *Blood Coagul Fibrinolysis.* 2014;25: 277–279.
- [57] Mugler K, Lefkowitz JB. False-positive D-dimer result in a patient with Castleman disease. *Arch Pathol Lab Med.* 2004;128:328–331.
- [58] Pittet JL, de Moerloose P, Reber G, et al. VIDAS D-dimer: fast quantitative ELISA for measuring D-dimer in plasma. *Clin Chem.* 1996;42:410–415.
- [59] La'ulu SL, Dominguez CM, Roberts WL. Performance characteristics of the AxsYM D-dimer assay. *Clin Chim Acta.* 2008;390:148–151.
- [60] Park SJ, Chi HS, Chun SH, et al. Evaluation of performance including influence by interfering substances of the innovance D-dimer assay on the Sysmex coagulation analyzer. *Clin Exp Thromb Hemost.* 2018; 4:11–14.
- [61] Chen L, Chen Y. Performance evaluation of the Sysmex CS-5100 automated coagulation analyzer. *Clin Lab.* 2015;61:653–660.
- [62] Flieder T, Gripp T, Knabbe C, et al. The Sysmex CS-5100 coagulation analyzer offers comparable analytical performance and excellent throughput capabilities. *Pract Lab Med.* 2016;6:38–47.
- [63] Huang H, Li H, Li D. Effect of serum monoclonal protein concentration on haemostasis in patients with multiple myeloma. *Blood Coagul Fibrinolysis.* 2015; 26:555–559.
- [64] Roller RE, Lahousen T, Lipp RW, et al. Elevated D-dimer results in a healthy patient. *Blood Coagul Fibrinolysis.* 2001;12:501–502.
- [65] Wu XY, Yin YF, Teng JL, et al. IgMk paraprotein from gammopathy patient can bind to cardiolipin and interfere with coagulation assay: a case report. *BMC Immunol.* 2017;18:32.
- [66] Gardiner C, Pennaneac'h C, Mackie IJ, et al. Falsely elevated D-dimer results in a healthy patient on account of heterophile antibodies. *Br J Haematol.* 2003;122:871–873.
- [67] Rouviere JA, Devignes J, de Maistre E, et al. Discrepancy between two methods of D-dimers measurement: one case of human anti-mouse antibody interference. *Ann Biol Clin (Paris).* 2008;66: 441–446.
- [68] Favresse J, Burlacu MC, Maiter D, et al. Interferences with thyroid function immunoassays: clinical

- implications and detection algorithm. *Endocr Rev.* 2018;39:830–850.
- [69] Bohm-Weigert M, Wissel T, Muth H, et al. Long- and short-term in vitro D-dimer stability measured with innovance D-Dimer. *Thromb Haemost.* 2010;103:461–465.
- [70] Betsou F, Roussel B, Guillaume N, et al. Long-term stability of coagulation variables: protein S as a biomarker for preanalytical storage-related variations in human plasma. *Thromb Haemost.* 2009;101:1172–1175.
- [71] Foshat M, Bates S, Russo W, et al. Effect of freezing plasma at -20 degrees C for 2 weeks on prothrombin time, activated partial thromboplastin time, dilute Russell viper venom time, activated protein C resistance, and D-dimer levels. *Clin Appl Thromb Hemost.* 2015;21:41–47.
- [72] Woodhams B, Girardot O, Blanco MJ, et al. Stability of coagulation proteins in frozen plasma. *Blood Coagul Fibrinolysis.* 2001;12:229–236.
- [73] Bastin P, Favresse J, Streel C, et al. Assessment of in vitro stability: a call for harmonization across studies. *Clin Chem Lab Med.* 2018;56:e121–e124.
- [74] Topic E, Nikolac N, Panteghini M, et al. How to assess the quality of your analytical method? *Clin Chem Lab Med.* 2015;53:1707–1718.
- [75] Zurcher M, Sulzer I, Barizzi G, et al. Stability of coagulation assays performed in plasma from citrated whole blood transported at ambient temperature. *Thromb Haemost.* 2008;99:416–426.
- [76] Kemkes-Matthes B, Fischer R, Peetz D. Influence of 8 and 24-h storage of whole blood at ambient temperature on prothrombin time, activated partial thromboplastin time, fibrinogen, thrombin time, antithrombin and D-dimer. *Blood Coagul Fibrinolysis.* 2011;22:215–220.
- [77] Rimac V, Coen Herak D. Is it acceptable to use coagulation plasma samples stored at room temperature and 4 degrees C for 24 hours for additional prothrombin time, activated partial thromboplastin time, fibrinogen, antithrombin, and D-dimer testing? *Int J Lab Hem.* 2017;39:475–481.
- [78] Zhao Y, Lv G. Influence of temperature and storage duration on measurement of activated partial thromboplastin time, D-dimers, fibrinogen, prothrombin time and thrombin time, in citrate-anticoagulated whole blood specimens. *Int Jnl Lab Hem.* 2013;35:566–570.
- [79] Caliezi C, Reber G, Lammle B, et al. Agreement of D-dimer results measured by a rapid ELISA (VIDAS) before and after storage during 24h or transportation of the original whole blood samples. *Thromb Haemost.* 2000;83:177–178.
- [80] Gosselin RC, Dwyre DW. Determining the effect of freezing on coagulation testing: comparison of results between fresh and once frozen-thawed plasma. *Blood Coagul Fibrinolysis.* 2015;26:69–74.
- [81] Rylatt DB, Blake AS, Cottis LE, et al. An immunoassay for human D dimer using monoclonal antibodies. *Thromb Res.* 1983;31:767–778.
- [82] Greenberg CS, Devine DV, McCrae KM. Measurement of plasma fibrin D-dimer levels with the use of a monoclonal antibody coupled to latex beads. *Am J Clin Pathol.* 1987;87:94–100.
- [83] Lippi G, Lima-Oliveira G, Brocco G, et al. Estimating the intra- and inter-individual imprecision of manual pipetting. *Clin Chem Lab Med.* 2017;55:962–966.
- [84] Lippi G, Ippolito L, Russello T, et al. Analytical performance of the new ACL AcuStar HemosIL D-Dimer. *Blood Coagul Fibrinolysis.* 2012;23:164–167.
- [85] Salvagno GL, Lippi G, Manzato F, et al. Analytical comparison of AxSYM, HemosIL DD HS and innovance D-dimer immunoassays with the Vidas D-dimer. *Int J Lab Hematol.* 2009;31:475–477.
- [86] van der Graaf F, van den Borne H, van der Kolk M, et al. Exclusion of deep venous thrombosis with D-dimer testing-comparison of 13 D-dimer methods in 99 outpatients suspected of deep venous thrombosis using venography as reference standard. *Thromb Haemost.* 2000;83:191–198.
- [87] Mountain D, Jacobs I, Haig A. The VIDAS D-dimer test for venous thromboembolism: a prospective surveillance study shows maintenance of sensitivity and specificity when used in normal clinical practice. *Am J Emerg Med.* 2007;25:464–471.
- [88] Perrier A, Desmarais S, Miron MJ, et al. Non-invasive diagnosis of venous thromboembolism in outpatients. *Lancet.* 1999;353:190–195.
- [89] Sidelmann JJ, Gram J, Larsen A, et al. Analytical and clinical validation of a new point-of-care testing system for determination of D-Dimer in human blood. *Thromb Res.* 2010;126:524–530.
- [90] Helmersson-Karlqvist J, Karlsson B, Fredriksson A, et al. Evaluation of the Alere D-dimer test for point of care testing. *J Thromb Thrombolysis.* 2014;38:250.
- [91] Lucassen WA, Erkens PM, Geersing GJ, et al. Qualitative point-of-care D-dimer testing compared with quantitative D-dimer testing in excluding pulmonary embolism in primary care. *J Thromb Haemost.* 2015;13:1004–1009.
- [92] Sen B, Kesteven P, Avery P. Comparison of D-dimer point of care test (POCT) against current laboratory test in patients with suspected venous thromboembolism (VTE) presenting to the emergency department (ED). *J Clin Pathol.* 2014;67:437–440.
- [93] Marquardt U, Apau D. Point-of-care D-dimer testing in emergency departments. *Emerg Nurse.* 2015;23:29–35.
- [94] Rooney KD, Schilling UM. Point-of-care testing in the overcrowded emergency department-can it make a difference? *Crit Care.* 2014;18:692.
- [95] Geersing GJ, Toll DB, Janssen KJ, et al. Diagnostic accuracy and user-friendliness of 5 point-of-care D-dimer tests for the exclusion of deep vein thrombosis. *Clin Chem.* 2010;56:1758–1766.
- [96] Dempfle CE, Zips S, Ergul H, et al. The fibrin assay comparison trial (FACT): correlation of soluble fibrin assays with D-dimer. *Thromb Haemost.* 2001;86:1204–1209.
- [97] Meijer P, Haverkate F, Kluft C, et al. A model for the harmonisation of test results of different quantitative D-dimer methods. *Thromb Haemost.* 2006;95:567–572.

- [98] Spannagl M, Haverkate F, Reinauer H, et al. The performance of quantitative D-dimer assays in laboratory routine. *Blood Coagul Fibrinolysis*. 2005;16: 439–443.
- [99] Olson JD, Cunningham MT, Higgins RA, et al. D-dimer: simple test, tough problems. *Arch Pathol Lab Med*. 2013;137:1030–1038.
- [100] Jennings I, Woods TA, Kitchen DP, et al. Laboratory D-dimer measurement: improved agreement between methods through calibration. *Thromb Haemost*. 2007;98:1127–1135.
- [101] Madoiwa S, Kitajima I, Ohmori T, et al. Distinct reactivity of the commercially available monoclonal antibodies of D-dimer and plasma FDP testing to the molecular variants of fibrin degradation products. *Thromb Res*. 2013;132:457–464.
- [102] Nieuwenhuizen W. A reference material for harmonisation of D-dimer assays. Fibrinogen Subcommittee of the Scientific and Standardization Committee of the International Society of Thrombosis and Haemostasis. *Thromb Haemost*. 1997;77:1031–1033.
- [103] Longstaff C, Adcock D, Olson JD, et al. Harmonisation of D-dimer – A call for action. *Thromb Res*. 2016;137: 219–220.
- [104] Dempfle CE. D-dimer assays: the current status and new assay technologies. *Thromb Res*. 2006;118: 569–571.
- [105] Lippi G, Tripodi A, Simundic AM, et al. International survey on D-dimer test reporting: a call for standardization. *Semin Thromb Hemost*. 2015;41:287–293.
- [106] Goodwin AJ, Higgins RA, Moser KA, et al. Issues surrounding age-adjusted d-dimer cutoffs that practicing physicians need to know when evaluating patients with suspected pulmonary embolism. *Ann Intern Med*. 2017;166:361–363.
- [107] Giannitsis E, Mair J, Christersson C, et al. How to use D-dimer in acute cardiovascular care. *Eur Heart J Acute Cardiovasc Care*. 2017;6:69–80.
- [108] Tripodi A, Chantarangkul V. Performance of quantitative D-dimer methods: results of the Italian external quality assessment scheme. *J Thromb Haemost*. 2007;5:184–185.
- [109] Kitchen S. Quality assurance in hemostasis. A European versus North American perspective. *Lab Hematol*. 2004;10:175–177.
- [110] Siragusa S. D-dimer testing: advantages and limitations in emergency medicine for managing acute venous thromboembolism. *Int Emergency Med*. 2006;1:59–66.
- [111] Keeling DM, Mackie IJ, Moody A, et al. The diagnosis of deep vein thrombosis in symptomatic outpatients and the potential for clinical assessment and D-dimer assays to reduce the need for diagnostic imaging. *Br J Haematol*. 2004;124:15–25.
- [112] Pernod G, Wu H, de Maistre E, et al. Validation of STA-Liatest D-Di assay for exclusion of pulmonary embolism according to the latest Clinical and Laboratory Standard Institute/Food and Drug Administration guideline. Results of a multicenter management study. *Blood Coagul Fibrinolysis*. 2017; 28:254–260.
- [113] Boeer K, Siegmund R, Schmidt D, et al. Comparison of six D-dimer assays for the detection of clinically suspected deep venous thrombosis of the lower extremities. *Blood Coagul Fibrinolysis*. 2009;20:141–145.
- [114] Coen Herak D, Milos M, Zadro R. Evaluation of the innovance D-DIMER analytical performance. *Clin Chem Lab Med*. 2009;47:945–951.
- [115] de Moerloose P, Palareti G, Aguilar C, et al. A multicenter evaluation of a new quantitative highly sensitive D-dimer assay for exclusion of venous thromboembolism. *Thromb Haemost*. 2008;100:505–512.
- [116] Lucassen W, Geersing GJ, Erkens PM, et al. Clinical decision rules for excluding pulmonary embolism: a meta-analysis. *Ann Intern Med*. 2011;155:448–460.
- [117] Deitelzweig SB, Johnson BH, Lin J, et al. Prevalence of clinical venous thromboembolism in the USA: current trends and future projections. *Am J Hematol*. 2011;86:217–220.
- [118] Schouten HJ, Geersing GJ, Koek HL, et al. Diagnostic accuracy of conventional or age adjusted D-dimer cut-off values in older patients with suspected venous thromboembolism: systematic review and meta-analysis. *BMJ*. 2013;346:f2492.
- [119] Douma RA, le Gal G, Sohne M, et al. Potential of an age adjusted D-dimer cut-off value to improve the exclusion of pulmonary embolism in older patients: a retrospective analysis of three large cohorts. *BMJ*. 2010;340:c1475.
- [120] Righini M, Van Es J, Den Exter PL, et al. Age-adjusted D-dimer cutoff levels to rule out pulmonary embolism: the ADJUST-PE study. *Jama*. 2014;311:1117–1124.
- [121] Pernod G, Maignan M, Marlu R. Questioning the use of an age-adjusted D-dimer threshold to exclude venous thromboembolism: analysis of individual patient data from two diagnostic studies: comment. *J Thromb Haemost*. 2016;14:2553–2554.
- [122] Monagle P, Barnes C, Ignjatovic V, et al. Developmental haemostasis. Impact for clinical haemostasis laboratories. *Thromb Haemost*. 2006;95: 362–372.
- [123] Szecsi PB, Jorgensen M, Klajnbard A, et al. Haemostatic reference intervals in pregnancy. *Thromb Haemost*. 2010;103:718–727.
- [124] Tromeur C, van der Pol LM, Klok FA, et al. Pitfalls in the diagnostic management of pulmonary embolism in pregnancy. *Thromb Res*. 2017;151(Suppl 1): S86–S91.
- [125] Hunt BJ, Parmar K, Horspool K, et al. The DIPEP (Diagnosis of PE in Pregnancy) biomarker study: An observational cohort study augmented with additional cases to determine the diagnostic utility of biomarkers for suspected venous thromboembolism during pregnancy and puerperium. *Br J Haematol*. 2018;180:694–704.
- [126] Lippi G, Bonfanti L, Saccenti C, et al. Causes of elevated D-dimer in patients admitted to a large urban emergency department. *Eur J Intern Med*. 2014;25: 45–48.
- [127] Tripodi A. D-dimer testing in laboratory practice. *Clin Chem*. 2011;57:1256–1262.
- [128] Garcia Suquia A, Alonso-Fernandez A, de la Pena M, et al. High D-dimer levels after stopping

- anticoagulants in pulmonary embolism with sleep apnoea. *Eur Respir J*. 2015;46:1691–1700.
- [129] Lippi G, Mattiuzzi C, Franchini M. Sleep apnea and venous thromboembolism. A systematic review. *Thromb Haemost*. 2015;114:958–963.
- [130] Saposnik G, Barinagarrementeria F, Brown RD, Jr, et al. Diagnosis and management of cerebral venous thrombosis: a statement for healthcare professionals from the American Heart Association/American Stroke Assoc. 2011;42:1158–1192.
- [131] Dentali F, Squizzato A, Marchesi C, et al. D-dimer testing in the diagnosis of cerebral vein thrombosis: a systematic review and a meta-analysis of the literature. *J Thromb Haemost*. 2012;10:582–589.
- [132] Alons IM, Jellema K, Wermer MJ, et al. D-dimer for the exclusion of cerebral venous thrombosis: a meta-analysis of low risk patients with isolated headache. *BMC Neurol*. 2015;15:118.
- [133] Cui JS, Jing ZP, Zhuang SJ, et al. D-dimer as a biomarker for acute aortic dissection: a systematic review and meta-analysis. *Medicine (Baltimore)*. 2015; 94:e471.
- [134] Hiratzka LF, Bakris GL, Beckman JA, et al. 2010 ACCF/AHA/AATS/ACR/ASA/SCA/SCAI/SIR/STS/SVM guidelines for the diagnosis and management of patients with thoracic aortic disease: a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines, American Association for Thoracic Surgery, American College of Radiology, American Stroke Association, Society of Cardiovascular Anesthesiologists, Society for Cardiovascular Angiography and Interventions, Society of Interventional Radiology, Society of Thoracic Surgeons, and Society for Vascular Medicine. *Circulation*. 2010;121:e266–e369.
- [135] Suzuki T, Distant A, Zizza A, et al. Diagnosis of acute aortic dissection by D-dimer: the International Registry of Acute Aortic Dissection Substudy on Biomarkers (IRAD-Bio) experience. *Circulation*. 2009; 119:2702–2707.
- [136] Shimony A, Fillon KB, Mottillo S, et al. Meta-analysis of usefulness of d-dimer to diagnose acute aortic dissection. *Am J Cardiol*. 2011;107:1227–1234.
- [137] Watanabe H, Horita N, Shibata Y, et al. Diagnostic test accuracy of D-dimer for acute aortic syndrome: systematic review and meta-analysis of 22 studies with 5000 subjects. *Sci Rep*. 2016;6:26893.
- [138] Nazerian P, Morello F, Vanni S, et al. Combined use of aortic dissection detection risk score and D-dimer in the diagnostic workup of suspected acute aortic dissection. *Int J Cardiol*. 2014;175:78–82.
- [139] Rogers AM, Hermann LK, Booher AM, et al. Sensitivity of the aortic dissection detection risk score, a novel guideline-based tool for identification of acute aortic dissection at initial presentation: results from the international registry of acute aortic dissection. *Circulation*. 2011;123:2213–2218.
- [140] Fang Z, Zhu XH, Wei X, et al. The diagnostic value of combined D-dimer with other indicators in suspected acute aortic dissection patients. *Int J Cardiol*. 2018;268:215.
- [141] Nazerian P, Mueller C, Soeiro AM, et al. Diagnostic accuracy of the aortic dissection detection risk score plus d-dimer for acute aortic syndromes: the ADVISED prospective multicenter study. *Circulation*. 2018;137:250–258.
- [142] Akyildiz H, Akcan A, Ozturk A, et al. The correlation of the D-dimer test and biphasic computed tomography with mesenteric computed tomography angiography in the diagnosis of acute mesenteric ischemia. *Am J Surg*. 2009;197:429–433.
- [143] Chiu YH, Huang MK, How CK, et al. D-dimer in patients with suspected acute mesenteric ischemia. *Am J Emerg Med*. 2009;27:975–979.
- [144] Sun Y, Wang L, Wei X, et al. Analysis of a Chinese pedigree with Zellweger syndrome reveals a novel PEX1 mutation by next-generation sequencing. *Clin Chim Acta*. 2013;417:57–61.
- [145] Yang S, Fan X, Ding W, et al. D-dimer as an early marker of severity in patients with acute superior mesenteric venous thrombosis. *Medicine (Baltimore)*. 2014;93:e270.
- [146] Hmoud B, Singal AK, Kamath PS. Mesenteric venous thrombosis. *J Clin Exp Hepatol*. 2014;4:257–263.
- [147] Tilsed JV, Casamassima A, Kurihara H, et al. ESTES guidelines: acute mesenteric ischaemia. *Eur J Trauma Emerg Surg*. 2016;42:253–270.
- [148] Yikilmaz A, Karahan OI, Senol S, et al. Value of multislice computed tomography in the diagnosis of acute mesenteric ischemia. *Eur J Radiol*. 2011;80: 297–302.
- [149] Sun DL, Li SM, Cen YY, et al. Accuracy of using serum D-dimer for diagnosis of acute intestinal ischemia: A meta-analysis. *Medicine (Baltimore)*. 2017;96: e6380.
- [150] Brotman DJ, Segal JB, Jani JT, et al. Limitations of D-dimer testing in unselected inpatients with suspected venous thromboembolism. *Am J Med*. 2003; 114:276–282.
- [151] Stein PD, Hull RD, Patel KC, et al. D-dimer for the exclusion of acute venous thrombosis and pulmonary embolism: a systematic review. *Ann Intern Med*. 2004;140:589–602.
- [152] Heit JA, Spencer FA, White RH. The epidemiology of venous thromboembolism. *J Thromb Thrombolysis*. 2016;41:3–14.
- [153] Mazzolai L, Aboyans V, Ageno W, et al. Diagnosis and management of acute deep vein thrombosis: a joint consensus document from the European society of cardiology working groups of aorta and peripheral circulation and pulmonary circulation and right ventricular function. *Eur Heart J*. 2017. DOI: 10.1093/eurheartj/ehx003.
- [154] Delluc A, Tromeur C, Le Ven F, et al. Current incidence of venous thromboembolism and comparison with 1998: a community-based study in Western France. *Thromb Haemost*. 2016;116:967–974.
- [155] Goldhaber SZ, Bounameaux H. Pulmonary embolism and deep vein thrombosis. *Lancet*. 2012;379: 1835–1846.
- [156] Heit JA, Silverstein MD, Mohr DN, et al. Predictors of survival after deep vein thrombosis and pulmonary

- embolism: a population-based, cohort study. *Arch Intern Med.* 1999;159:445–453.
- [157] White RH. The epidemiology of venous thromboembolism. *Circulation.* 2003;107:14–18.
- [158] Spyropoulos AC, Lin J. Direct medical costs of venous thromboembolism and subsequent hospital readmission rates: an administrative claims analysis from 30 managed care organizations. *Jmcp.* 2007;13:475–486.
- [159] Jha AK, Larizgoitia I, Audera-Lopez C, et al. The global burden of unsafe medical care: analytic modelling of observational studies. *BMJ Qual Saf.* 2013;22:809–815.
- [160] Cohoon KP, Leibson CL, Ransom JE, et al. Direct medical costs attributable to venous thromboembolism among persons hospitalized for major operation: a population-based longitudinal study. *Surgery.* 2015;157:423–431.
- [161] Cohoon KP, Leibson CL, Ransom JE, et al. Costs of venous thromboembolism associated with hospitalization for medical illness. *Am J Manag Care.* 2015;21:e255–e263.
- [162] Grosse SD. Incidence-based cost estimates require population-based incidence data. A critique of Mahan et al. *Thromb Haemost.* 2011;107:192–193.
- [163] CLSI. Quantitative D-dimer for exclusion of venous thromboembolic disease. CLSI document H59-A. Wayne, PA: Clinical and Laboratory Standard Institute; 2011.
- [164] Weitz JJ, Fredenburgh JC, Eikelboom JW. A test in context: D-dimer. *J Am Coll Cardiol.* 2017;70:2411–2420.
- [165] Di Nisio M, Squizzato A, Rutjes AW, et al. Diagnostic accuracy of D-dimer test for exclusion of venous thromboembolism: a systematic review. *J Thromb Haemost.* 2007;5:296–304.
- [166] Roy PM, Colombet I, Durieux P, et al. Systematic review and meta-analysis of strategies for the diagnosis of suspected pulmonary embolism. *BMJ.* 2005;331:259.
- [167] Carrier M, Righini M, Djurabi RK, et al. VIDAS D-dimer in combination with clinical pre-test probability to rule out pulmonary embolism. A systematic review of management outcome studies. *Thromb Haemost.* 2009;101:886–892.
- [168] Mullier F, Vanpee D, Jamart J, et al. Comparison of five D-dimer reagents and application of an age-adjusted cut-off for the diagnosis of venous thromboembolism in emergency department. *Blood Coagul Fibrinolysis.* 2014;25:309–315.
- [169] Geersing GJ, Janssen KJ, Oudega R, et al. Excluding venous thromboembolism using point of care D-dimer tests in outpatients: a diagnostic meta-analysis. *BMJ.* 2009;339:b2990.
- [170] Chunilal SD, Eikelboom JW, Attia J, et al. Does this patient have pulmonary embolism? *Jama.* 2003;290:2849–2858.
- [171] Jaeschke R, Guyatt GH, Sackett DL. Users' guides to the medical literature. III. How to use an article about a diagnostic test. B. What are the results and will they help me in caring for my patients? The Evidence-Based Medicine Working Group. *Jama.* 1994;271:703–707.
- [172] Fagan TJ. Letter: nomogram for Bayes theorem. *N Engl J Med.* 1975;293:257.
- [173] Penalzoza A, Roy PM, Kline J, et al. Performance of age-adjusted D-dimer cut-off to rule out pulmonary embolism. *J Thromb Haemost.* 2012;10:1291–1296.
- [174] Righini M, Perrier A, De Moerloose P, et al. D-Dimer for venous thromboembolism diagnosis: 20 years later. *J Thromb Haemost.* 2008;6:1059–1071.
- [175] Gardiner C, Pennaneac'h C, Walford C, et al. An evaluation of rapid D-dimer assays for the exclusion of deep vein thrombosis. *Br J Haematol.* 2005;128:842–848.
- [176] Wells PS, Hirsh J, Anderson DR, et al. A simple clinical model for the diagnosis of deep-vein thrombosis combined with impedance plethysmography: potential for an improvement in the diagnostic process. *J Intern Med.* 1998;243:15–23.
- [177] Robert-Ebadi H, Mostaguir K, Hovens MM, et al. Assessing clinical probability of pulmonary embolism: prospective validation of the simplified Geneva score. *J Thromb Haemost.* 2017;15:1764–1769.
- [178] Spencer FA, Goldberg RJ, Lessard D, et al. Factors associated with adverse outcomes in outpatients presenting with pulmonary embolism: the Worcester Venous Thromboembolism Study. *Circ Cardiovasc Qual Outcomes.* 2010;3:390.
- [179] Penalzoza A, Verschuren F, Meyer G, et al. Comparison of the unstructured clinician gestalt, the wells score, and the revised Geneva score to estimate pretest probability for suspected pulmonary embolism. *Ann Emerg Med.* 2013;62:117–124 e2.
- [180] van der Hulle T, Cheung WY, Kooij S, et al. Simplified diagnostic management of suspected pulmonary embolism (the years study): a prospective, multi-centre, cohort study. *Lancet.* 2017;390:289–297.
- [181] Kristoffersen AH, Ajzner E, Bauca JM, et al. Pre- and post-test probabilities of venous thromboembolism and diagnostic accuracy of D-dimer, estimated by European clinicians working in emergency departments. *Thromb Res.* 2017;159:19–23.
- [182] Kristoffersen AH, Ajzner E, Rogic D, et al. Is D-dimer used according to clinical algorithms in the diagnostic work-up of patients with suspicion of venous thromboembolism? A study in six European countries. *Thromb Res.* 2016;142:1–7.
- [183] Le Gal G, Righini M, Wells PS. D-dimer for pulmonary embolism. *Jama.* 2015;313:1668–1669.
- [184] Konstantinides SV, Torbicki A, Agnelli G, et al. 2014 ESC guidelines on the diagnosis and management of acute pulmonary embolism. *Eur Heart J.* 2014;35:3033–3069.
- [185] Linkins LA, Bates SM, Lang E, et al. Selective D-dimer testing for diagnosis of a first suspected episode of deep venous thrombosis. *Ann Intern Med.* 2013;158:93–100.
- [186] Pasha SM, Klok FA, Snoep JD, et al. Safety of excluding acute pulmonary embolism based on an unlikely clinical probability by the Wells rule and normal D-dimer concentration: a meta-analysis. *Thromb Res.* 2010;125:e123–e127.
- [187] Roy PM, Meyer G, Vielle B, et al. Appropriateness of diagnostic management and outcomes of suspected

- pulmonary embolism. *Ann Intern Med.* 2006;144:157–164.
- [188] Bates SM, Takach Lapner S, Douketis JD, et al. Rapid quantitative D-dimer to exclude pulmonary embolism: a prospective cohort management study. *J Thromb Haemost.* 2016;14:504–509.
- [189] Bates SM. D-dimer assays in diagnosis and management of thrombotic and bleeding disorders. *Semin Thromb Hemost.* 2012;38:673–682.
- [190] Sohne M, Kamphuisen PW, van Mierlo PJ, et al. Diagnostic strategy using a modified clinical decision rule and D-dimer test to rule out pulmonary embolism in elderly in- and outpatients. *Thromb Haemost.* 2005;94:206–210.
- [191] Farm M, Siddiqui A, Onelöv L, et al. Age-adjusted D-dimer cutoffs for DVT and pulmonary embolism: a comparison of five assays. *Blood.* 2016;128:1419–1419.
- [192] Prochaska JH, Frank B, Nagler M, et al. Age-related diagnostic value of D-dimer testing and the role of inflammation in patients with suspected deep vein thrombosis. *Sci Rep.* 2017;7:4591.
- [193] Raja AS, Greenberg JO, Qaseem A, et al. Evaluation of patients with suspected acute pulmonary embolism: best practice advice from the Clinical Guidelines Committee of the American College of Physicians. *Ann Intern Med.* 2015;163:701–711.
- [194] Solberg R, Glass G. Adjusting D-dimer cutoffs: brief literature summary and issues in clinical use. *Am J Emerg Med.* 2018;36:2105–2107.
- [195] Kabrhel C, Mark Courtney D, Camargo CA, Jr, et al. Potential impact of adjusting the threshold of the quantitative D-dimer based on pretest probability of acute pulmonary embolism. *Acad Emerg Med.* 2009;16:325–332.
- [196] Linkins LA, Bates SM, Ginsberg JS, et al. Use of different D-dimer levels to exclude venous thromboembolism depending on clinical pretest probability. *J Thromb Haemost.* 2004;2:1256–1260.
- [197] Takach Lapner S, Julian JA, Linkins LA, et al. Comparison of clinical probability-adjusted D-dimer and age-adjusted D-dimer interpretation to exclude venous thromboembolism. *Thromb Haemost.* 2017;117:1937–1943.
- [198] Kline JA, Hogg MM, Courtney DM, et al. D-dimer threshold increase with pretest probability unlikely for pulmonary embolism to decrease unnecessary computerized tomographic pulmonary angiography. *J Thromb Haemost.* 2012;10:572–581.
- [199] van der Pol LM, van der Hulle T, Mairuhu ATA, et al. Combination of pulmonary embolism rule-out criteria and YEARS algorithm in a European cohort of patients with suspected pulmonary embolism. *Thromb Haemost.* 2018;118:547–552.
- [200] Desciak MC, Martin DE. Perioperative pulmonary embolism: diagnosis and anesthetic management. *J Clin Anesth.* 2011;23:153–165.
- [201] Penalzoa AKC, Verschuren F, Vielle B, et al. D-dimer relevance score. *J Thromb Haemost.* 2013;11:1–1322.
- [202] Wang M, Lu S, Li S, et al. Reference intervals of D-dimer during the pregnancy and puerperium period on the STA-R evolution coagulation analyzer. *Clin Chim Acta.* 2013;425:176–180.
- [203] Kline JA, Williams GW, Hernandez-Nino J. D-dimer concentrations in normal pregnancy: new diagnostic thresholds are needed. *Clin Chem.* 2005;51:825–829.
- [204] Kovac M, Mikovic Z, Rakicevic L, et al. The use of D-dimer with new cutoff can be useful in diagnosis of venous thromboembolism in pregnancy. *Eur J Obstet Gynecol Reprod Biol.* 2010;148:27–30.
- [205] Van der Pol LM, Mairuhu AT, Tromeur C, et al. Use of clinical prediction rules and D-dimer tests in the diagnostic management of pregnant patients with suspected acute pulmonary embolism. *Blood Rev.* 2017;31:31–36.
- [206] Leung AN, Bull TM, Jaeschke R, et al. American Thoracic Society documents: an official American Thoracic Society/Society of Thoracic Radiology Clinical Practice Guideline-Evaluation of Suspected Pulmonary Embolism in Pregnancy. *Radiology.* 2012;262:635–646.
- [207] Sökücü S, Gökçe S, Güllüoğlu M, et al. The role of the non-invasive serum marker FibroTest-ActiTest in the prediction of histological stage of fibrosis and activity in children with naïve chronic hepatitis B infection. *Scand J Infect Dis.* 2010;42:699–703.
- [208] Geersing GJ, Zuithoff NP, Kearon C, et al. Exclusion of deep vein thrombosis using the Wells rule in clinically important subgroups: individual patient data meta-analysis. *BMJ.* 2014;348:g1340.
- [209] Carrier M, Lee AY, Bates SM, et al. Accuracy and usefulness of a clinical prediction rule and D-dimer testing in excluding deep vein thrombosis in cancer patients. *Thromb Res.* 2008;123:177–183.
- [210] Ay C, Dunkler D, Marosi C, et al. Prediction of venous thromboembolism in cancer patients. *Blood.* 2010;116:5377–5382.
- [211] Robert-Ebadi H, Bertoletti L, Combescure C, et al. Effects of impaired renal function on levels and performance of D-dimer in patients with suspected pulmonary embolism. *Thromb Haemost.* 2014;112:614–620.
- [212] Lindner G, Funk GC, Pfortmueller CA, et al. D-dimer to rule out pulmonary embolism in renal insufficiency. *Am J Med.* 2014;127:343–347.
- [213] Karami-Djurabi R, Klok FA, Kooiman J, et al. D-dimer testing in patients with suspected pulmonary embolism and impaired renal function. *Am J Med.* 2009;122:1050–1053.
- [214] Douketis J, Toseito A, Marcucci M, et al. Risk of recurrence after venous thromboembolism in men and women: patient level meta-analysis. *BMJ.* 2011;342:d813.
- [215] Palareti G, Cosmi B, Legnani C, et al. D-dimer testing to determine the duration of anticoagulation therapy. *N Engl J Med.* 2006;355:1780–1789.
- [216] Legnani C, Palareti G, Cosmi B, et al. Different cut-off values of quantitative D-dimer methods to predict the risk of venous thromboembolism recurrence: a post-hoc analysis of the PROLONG study. *Haematologica.* 2008;93:900–907.
- [217] Legnani C, Cini M, Cosmi B, et al. Age and gender specific cut-off values to improve the performance of

- D-dimer assays to predict the risk of venous thromboembolism recurrence. *Intern Emerg Med*. 2013;8:229–236.
- [218] Lippl G, Favalaro EJ, Cervellin G. A review of the value of D-dimer testing for prediction of recurrent venous thromboembolism with increasing age. *Semin Thromb Hemost*. 2014;40:634–639.
- [219] Verhovsek M, Douketis JD, Yi Q, et al. Systematic review: D-dimer to predict recurrent disease after stopping anticoagulant therapy for unprovoked venous thromboembolism. *Ann Intern Med*. 2008;149:481–490.
- [220] Nagler M, Ten Cate H, Prins MH, et al. Risk factors for recurrence in deep vein thrombosis patients following a tailored anticoagulant treatment incorporating residual vein obstruction. *Res Pract Thromb Haemost*. 2018;2:299–309.
- [221] Baglin T, Luddington R, Brown K, et al. Incidence of recurrent venous thromboembolism in relation to clinical and thrombophilic risk factors: prospective cohort study. *Lancet*. 2003;362:523–526.
- [222] Tosetto A, Iorio A, Marcucci M, et al. Predicting disease recurrence in patients with previous unprovoked venous thromboembolism: a proposed prediction score (DASH). *J Thromb Haemost*. 2012;10:1019–1025.
- [223] Tosetto A, Testa S, Martinelli I, et al. External validation of the DASH prediction rule: a retrospective cohort study. *J Thromb Haemost*. 2017;15:1963–1970.
- [224] Ensor J, Riley RD, Moore D, et al. Systematic review of prognostic models for recurrent venous thromboembolism (VTE) post-treatment of first unprovoked VTE. *BMJ Open*. 2016;6:e011190.
- [225] Eichinger S, Heinze G, Jandeck LM, et al. Risk assessment of recurrence in patients with unprovoked deep vein thrombosis or pulmonary embolism: the Vienna prediction model. *Circulation*. 2010;121:1630–1636.
- [226] Marcucci M, Iorio A, Douketis JD, et al. Risk of recurrence after a first unprovoked venous thromboembolism: external validation of the Vienna Prediction Model with pooled individual patient data. *J Thromb Haemost*. 2015;13:775–781.
- [227] Rodger MA, Le Gal G, Anderson DR, et al. Validating the HERDOO2 rule to guide treatment duration for women with unprovoked venous thrombosis: multinational prospective cohort management study. *BMJ*. 2017;356:j1065.
- [228] Rodger MA, Kahn SR, Wells PS, et al. Identifying unprovoked thromboembolism patients at low risk for recurrence who can discontinue anticoagulant therapy. *Cmaj*. 2008;179:417–426.
- [229] Ageno W, Squizzato A, Wells PS, et al. The diagnosis of symptomatic recurrent pulmonary embolism and deep vein thrombosis: guidance from the SSC of the ISTH. *J Thromb Haemost*. 2013;11:1597–1602.
- [230] Gibson CM, Jennings LK, Chi G, et al. Association of D-dimer levels with clinical event rates and the efficacy of betrixaban versus enoxaparin in the APEX trial. *TH Open*. 2018;02:e16–e24.
- [231] Gibson C, Spyropoulos A, Cohen A, et al. The IMPROVEDD VTE risk score: Incorporation of D-Dimer into the IMPROVE score to improve venous thromboembolism risk stratification. *TH Open*. 2017;01:e56–e65.
- [232] Cohen AT, Harrington R, Goldhaber SZ, et al. The design and rationale for the acute medically ill venous thromboembolism prevention with extended duration betrixaban (APEX) study. *Am Heart J*. 2014;167:335–341.
- [233] Desjardins L, Bara L, Boutitie F, et al. Correlation of plasma coagulation parameters with thromboprophylaxis, patient characteristics, and outcome in the MEDENOX study. *Arch Pathol Lab Med*. 2004;128:519–526.
- [234] Fan J, Li X, Cheng Y, et al. Measurement of D-dimer as aid in risk evaluation of VTE in elderly patients hospitalized for acute illness: a prospective, multicenter study in China. *Cim*. 2011;34:96–104.
- [235] Wada H, Thachil J, Di Nisio M, et al. Guidance for diagnosis and treatment of DIC from harmonization of the recommendations from three guidelines. *J Thromb Haemost*. 2013;11:761–767.
- [236] Wada H, Takahashi H, Uchiyama T, et al. The approval of revised diagnostic criteria for DIC from the Japanese Society on Thrombosis and Hemostasis. *Thromb J*. 2017;15:17.
- [237] Hayakawa M, Gando S, Hoshino H. A Prospective comparison of new Japanese criteria for disseminated intravascular coagulation: new Japanese criteria versus ISTH criteria. *Clin Appl Thromb Hemost*. 2007;13:172–181.
- [238] Asakura H, Takahashi H, Uchiyama T, et al. Proposal for new diagnostic criteria for DIC from the Japanese Society on Thrombosis and Hemostasis. *Thromb J*. 2016;14:42.
- [239] Levi M, Toh CH, Thachil J, et al. Guidelines for the diagnosis and management of disseminated intravascular coagulation. *British Committee for Standards in Haematology. Br J Haematol*. 2009;145:24–33.
- [240] Gris JC, Bouvier S, Cochery-Nouvellon E, et al. Fibrin-related markers in patients with septic shock: individual comparison of D-dimers and fibrin monomers impacts on prognosis. *Thromb Haemost*. 2011;106:1228–1230.

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Chapter 34

D-Dimer and Fibrinogen/Fibrin Degradation Products

Linda J. Stang

Abstract

Although clinical requests for D-dimer are generally in the minority of assays in the routine clinical laboratory, they are an important aspect—especially if the laboratory supports an active emergency room and hematology service. Throughout the literature, D-dimer assays have been used for many purposes in the research setting; however it is generally the negative predictive value of the assay that is the most common piece of information being utilized from the standpoint of a clinician. Research or clinical needs will dictate the type of assay required—a qualitative, semiquantitative, or quantitative D-dimer assay may be appropriate for a particular purpose. Commonalities and differences between these assay types are outlined here, as well as universal concerns regarding standardization of D-dimer assay results.

Key words D-dimer, Fibrinogen/fibrin degradation products, Quantitative, Semiquantitative, Qualitative, Standardization

1 Introduction

D-dimer assays are used in the clinical setting mainly to aid in diagnosis of deep vein thrombosis (DVT) and pulmonary embolus (PE) (see Notes 1 and 2). These assays are also useful for detection of post-op thrombosis, and quantitative or semiquantitative assays (to a lesser degree) can be used as a continuous marker of thrombolytic activity in a patient who has been diagnosed with arterial and/or venous thrombosis. Elevated levels can also be seen in, but are not limited to, disseminated intravascular coagulation (DIC) and trauma, as well as liver disease, neoplasia, inflammation, or even normal pregnancy (1, 2). Increased levels have also, more recently, been associated with increased risk of a recurrent thrombosis following anticoagulation therapy (2).

In order for a patient's sample to have elevated levels of D-dimer present, a number of reactions must have taken place *in vivo* prior to sample collection. It all begins simplistically with fibrinogen being acted upon by thrombin, effectively releasing fibrinopeptide A (FPA),


and then fibrinopeptide B to produce fibrin monomer, which progresses to a polymer. Activated factor XIII can then begin to stabilize the polymer, by exchanging the hydrogen bonds for covalent cross-links (3). This stabilized fibrin network is almost immediately attacked for subsequent degradation. Plasminogen that has been activated to plasmin, via its primary activator—tissue plasminogen activator (tPA) (4), has the ability to cleave both the fibrin contained in the clot and its circulating precursor, fibrinogen. The fibrin mesh can bind both the tPA and plasminogen, and when bound, the rate of activation increases substantially. Through this mechanism, the plasmin is located right where it is needed—at the site of the fibrin clot, as opposed to being free in the circulating plasma (5). Only the fibrin that has completed the cross-linking step can produce fragments with the actual D-dimer configuration. When broken down, native fibrinogen produces fragments, although none of them share this unique D-dimer identity. Therefore, the fragments produced from fibrinogen alone should not cross-react with assay reagents containing antibodies that are specific for the D-dimer neo-epitope. The cross-linked polymer is of vast size, at least on a molecular level, and the plasmin can attack anywhere that D and E domains occur together (within the original fibrinogen molecule configuration), producing a heterogeneity of fragments. Of the cross-linked degradation products formed, the smallest fragment type that is resistant to further plasmin degradation is actually denoted as a DDE fragment and this fragment size alone is sometimes referred to as D-dimer, when in actuality all fragments containing the dimeric D domain collectively are D-dimer fragments.

1.1 Nomenclature and Terminology for Fibrinogen, Fibrin, and Fragments

The following is a simplified version of nomenclature and terminology commonly used in literature and texts to discuss the process of fibrin production and the fragments produced (4, 6).

A native fibrinogen molecule is usually referred to as (D-E-D); within this molecule the D domains, which are at opposite ends of the molecule, are not dimeric (or adjacent to one another).

The action of thrombin and F XIIIa on this molecule produces a stable polymer denoted as

Plasmin has the capability to then act on this polymer, causing a breakage to occur only at the junctions denoted as (-) between the (D-E). Random breakage creates the possibility of the following fragments being generated, as well as various larger fragments:

Fragment Composition:	DD-E E-DD	E-DD-E DD-E-DD	DD E
Fragment name:	YD/DY	YY/DXD	DD/E

All of the degradation products are collectively referred to as D-dimer. Realistically, the lysis of a fibrin clot produces a very heterogeneous group of degradation products. These include X oligomers, YY/DXD, YD/DY, and DD/E, where (4):

Fragment X=D-E-D; an intermediate degradation product consisting of all three domain regions (resembles fibrinogen, but lacks fibrinopeptide A and B).

Fragment Y=D-E; contains a central E and one of either terminal D domains.

Initially large degradation products are released and only after a prolonged incubation with plasmin would the majority of fragments be DD/E, which is the smallest fragment resistant to further breakdown by plasmin (4). This is especially important when you consider the following: if a random group of patients were being tested for D-dimer, they would all be at different stages in the process of lysis. This fact, along with the different epitope specificity for the various antibodies in commercial kits, can produce slight variability in results between different assays.

Some publications also use the term fibrin degradation products (FDP), which occasionally can denote fibrinogen degradation products, as well. Another term sometimes used is XL-FDP, assumably in reference to "cross-linked FDP."

1.2 Detection of Fibrin/Fibrinogen Degradation Products

Fragments caused from plasmin breakdown of fibrinogen and fibrin can be detected using a commercially available fibrin/fibrinogen degradation product kit (Thrombo-Wellcotest, Remel, UK). This assay uses antibodies to the individual domains of D and E. Therefore the sample must be free of fibrinogen which *will* cross-react. A special collection tube containing *Bothrops atrox* venom and soya bean trypsin inhibitor will promote complete clotting of the fibrinogen, as well no further increase of D-dimer fragments should occur from the fibrin in the collection tube. Although not commonly used, this assay allows for testing for the breakdown products of fibrinogen and fibrin and has been particularly useful in DIC (7). When no pathological clot is present in other patient populations, this assay may be useful for detecting hyperactivity in a patient's fibrinolysis pathway. More specific assays relating to the root cause of the hyperactivity are available in specialized reference laboratories (plasminogen, tPA, PAI-1, etc.) and may be useful in these situations.

1.3 Differences in Assay Design

In order to choose an appropriate assay, the size of laboratory, type of equipment needed, acuity of patient population, and estimated frequency of the demand for D-dimer testing must be taken into consideration. Assays can be quantitative, semiquantitative, or qualitative (see Note 3). Quantitative assays would be more useful in situations where risk stratification is important or for ongoing monitoring in conditions such as DIC. The overall size of the measuring range of the assay is also important for this patient group as

patient results can reach substantial values. For the primary purpose of exclusion of PE or DVT, a qualitative assay or one with a small measuring range is acceptable (8).

Commercially available D-dimer assays can be broken down roughly into four categories, all of which rely on antibodies and the ensuing agglutination, to assist in detection:

1. Slide tests:

- (a) Qualitative slide tests, usually known for having a very high negative predictive value; suitable for emergency room triage or similar situations where a negative or a positive result is adequate. Assay principles are somewhat similar to standard enzyme-linked immunosorbent assays (ELISA), using a modified sandwich antibody technique, or a modified antibody with essentially two Fab portions which allow the reaction to occur.
- (b) Semiquantitative assays for which general ranges of D-dimer are given as results (i.e., <0.5 , $0.5-2.0$, $2.0-4.0$, and >4.0 $\mu\text{g/mL}$). The assay usually has a specific threshold over which a positive result will occur with a neat sample. The semiquantitative result is calculated by multiplying this threshold value by the reciprocal of the highest dilution of the sample which still produces a positive result. Slight modifications of the dilutions can usually be made to accommodate desired result ranges for these assays. When multiple dilutions are recommended by the manufacturer for more defined semiquantitative ranges, it may be more economic to omit one or more of the "middle" dilutions and only test the lower and higher of the recommended dilutions when a patient's neat sample is positive. Possibly at least one dilution can be omitted from the protocol without having any clinically significant differences produced in the semiquantitative result.

Modifications to a qualitative slide test are sometimes included in the manufacturer's instructions to enable semiquantitative analysis using the same kit. Usually the instructions involve running preset-specific dilutions of patient's plasma to enable determination of the range in which a patient's D-dimer value falls.

2. Automated quantitative methods: These are usually immunoturbidimetric assays. Most automated methods are designed for the mid- to large-size coagulation analyzers, utilizing a monoclonal (or a combination of two monoclonals) antibody-coated (usually latex) microbead suspension, which aggregates proportionally to the concentration of D-dimer present in the patient's sample, producing an increase in light scatter or a decrease in transmitted light. One advantage to these assays is that quantitative results are available in minutes.

3. Quantitative ELISA: This assay type accommodates large batch analysis and is suitable for research studies (some commercial products are labelled "For research use only"). Standard ELISA methodology usually takes a few hours to generate results, making this assay type unsuitable in some clinical situations, where assay results are requested on a stat basis (see Note 4).
4. Small device or point-of-care quantitative analyzers: Suitable for smaller sites that desire a quantitative result; usually these sites use this result as a triage tool, transferring patients with positive results to larger institutions for active treatment.

1.4 Sample Requirements

Most standard assays for D-dimer contain antibodies that do not react with native fibrinogen; therefore plasma is a suitable specimen type. The commercially available assay that does contain antibodies which cross-react with fibrinogen requires a special serum collection tube, which is discussed in Subheading 1.2.

Platelet-poor plasma separated from a 3.2% sodium citrate vacuum collection tube (centrifuged at $1,700 \times g$ for 10 min at room temperature) is generally the specimen required for a D-dimer assay. There are some variations to this, for particular commercially available assays. Some point-of-care analyzers require and some ELISA methods can substitute an EDTA (mauve top) sample in place of the citrate sample. The stability of an EDTA sample may be shorter, as some nonspecific agglutination may occur 4 h post collection (9). For some commercial assays, heparin primary collection tubes can also be a suitable sample replacement. It is important to verify and comply with the product insert specimen requirements for each specific assay.

Any Na citrate, EDTA, or heparin samples that are clotted should be discarded and recollected, as *ex vivo* sample clotting has the potential to cause a false positive result.

The sample stability for various commercially available D-dimer assays ranges from 4 to 8 h at room temperature, 2 weeks at -20°C , and 6 months at -70°C . Frozen samples should be vortexed after being thawed in a 37°C water bath for approximately 5–10 min, the length of thaw time being dependent on sample volume.

The *in vivo* half-life of D-dimer is approximately 12 h; thus timing of sample collection is not a critical issue, but this may be useful information when interpreting sequential quantitative or semi-quantitative results on the same patient over a given time span.

1.5 Difference in Antibody Specificity

Various kits use different clones of monoclonal antibodies, or combinations of clones of monoclonal antibodies. As a general rule, these antibodies have no cross-reactivity to the precursor molecule fibrinogen, FPA, or fibrin monomer.

The degradation products from fibrin clot lysis are a heterogeneous group of substances, as discussed earlier. Usually the term "D-dimer" is a collective term, in reference to this entire group of

fragments. D-dimer-E, YD, DXD, and X oligomers are the individual items that make up this group. Specificity of the various antibodies used in commercial assays does cause slight variation between assays. Usually these variations are limited to the region of low positive values, or those near the predetermined cutoffs for clinical investigation. For the most part any comparison of patient samples containing a more robust amount of D-dimer will give results corresponding in the appropriate range (relative to the various reference ranges for those kits) on comparison between various commercial assays. As there currently is no international standard for D-dimer, each kit and therefore each antibody (or antibody combination) must be clinically validated and the values determined for a particular assay *cannot* be interchanged with another, without appropriate comparison studies.

Usually manufacturers have done robust studies to determine appropriate cutoffs associated with clinical situations (i.e., negative predictive value for DVT). It is a time-consuming and overwhelming project to replicate these studies and is not usually feasible in most clinical laboratories. Small in-house studies, that would alter the predetermined cutoff used for commercially available products, are not recommended, as they usually do not have the mathematical power necessary to convey significance (10).

1.6 Standardization of Assays, Including Fibrinogen Equivalent Units

Fibrinogen equivalent units or FEU are not referred to in the majority of commercially available package inserts. This nomenclature is somewhat historical in nature. When D-dimer was first quantitated in various laboratories, the vast differences in methodology (including Staphylococcus clumping assay or tanned red cell hemagglutination inhibition assay (6)) produced a large difference in reference ranges and clinically significant values amongst the various assays. In an effort to standardize the assays, the FEU was implemented. One FEU is the quantity of fibrinogen initially present that leads to the observed level of D-dimer. As the D-dimer fragments are a heterogeneous group, and each of the historical assays relied on a completely different principle, the FEU helped to make assays comparable. For the purposes of practicality, a positive cutoff of 0.5 $\mu\text{g/mL}$ for FEU is generally equivalent to 0.25 $\mu\text{g/mL}$ of D-dimer (11).

One of the problems plaguing this assay is that there is currently no reference standard available for antigenic D-dimer (12). Various manufacturers of commercially available products also cite this lack of an international reference standard for D-dimer contributing to the usage of FEU in some assay kits. In the determination of D-dimer standards by various manufacturers, the degree to which plasmin is allowed to act on the preset amount of fibrinogen initially measured for the standard may also play a role in the difference between commercially available products. Some manufacturers also use pools of affected patients, for standardization of the assay. High-molecular-weight fractions may also be more indicative of fibrinolysis occurring in a zone of flow (13) and can be associated

with patients in DIC (14). To add to this already complex standardization situation, various manufacturers also do not use one common measurement unit for reporting. Some multi-assay comparison and standardization trials have increased our awareness and understanding of this complex problem (15); however at the time of this writing, the issue is still outstanding for the most part and must be considered when comparing data from various commercially available products.

2 Materials

The following is a table (Table 1) of common commercially available products for D-dimer and FDP. All required materials come with each kit. Patient values that fall below cutoff values are considered normal; for patients with levels above, suitable follow-up testing is usually recommended by the manufacturer (see Notes 1–6). Although the author attempted to include kits from all companies, this table may not be all inclusive for every product on the market. The intended objective is to leave an impression of the variability between assay values for commercial kits.

3 Methods

3.1 Slide Tests: Method and Procedural Notes

One of the downfalls of slide tests is the variability in the interpretation of the presence of agglutination by the end user; this is especially prominent if the lab does not do a large volume of slide tests. It is important to have good lighting by which to read the test and adhere to the product insert mixing instructions exactly. Some users have found that interpreting the agglutination present is easier if read from beneath a clear slide when the reaction is occurring; however this does present a safety hazard. Adequate proficiency testing for all staff performing the slide tests should be implemented to ensure quality results.

3.2 Method for Qualitative Slide Test

Qualitative slide methodologies vary between manufacturer with regard to sample amounts, speed and duration of slide rocking etc. Always use manufacturer's recommendations. For basic qualitative slide test instructions (see Table 2).

3.3 Semiquantitative Modification of Qualitative Slide Tests

The majority of qualitative slide tests can be modified to produce semiquantitative results. Serial (or other) dilutions of patient sample are made and tested in the same manner as for the qualitative slide test method above. The value of the patient's D-dimer is calculated using the threshold value for a positive result (on a neat patient sample), multiplied by the reciprocal of the dilution for the highest dilution of a patient that still tests positive in the reaction. For instance if the threshold value is 0.500 µg/mL and the highest

Table 1
Reference ranges and alert values

Product	Reference range or upper end of normal (NOTE: unit values differ)	Cutoff value for exclusion of DVT and/or PE, if different than reference range value
<i>Qualitative or semiquantitative assays</i>		
Clearview Simplify D-Dimer	<80 ng/mL	>80 ng/mL
AGEN Simplired D-Dimer	Not specified	>0.12 mg/mL
Diagnostica Stago D-Di Test	<0.5 µg/mL (FEU)	Not specified
HemosIL Dimertest	<0.20 mg/mL	Not specified
AMAX Accuclot D-Dimer	8–135 ng/mL	Not specified
Trinity Biotech Minutex D-Dimer	8–135 ng/mL	Not specified
Thrombo-Wellcotest FDP (serum sample)	4.9 ± 2.8 µg/mL	DVT—10–40 µg/mL PE—>40 µg/mL
Remel D-Dimer	<1 µg/mL	Not specified
<i>Automated quantitative assays</i>		
HemosIL D-Dimer HS	243 ng/mL	230 ng/mL
HemosIL D-Dimer (quant)	Dependent on analyzer models	230 ng/mL
Diagnostica Stago Liatest D-Di	<0.500 µg/mL (FEU)	>0.500 µg/mL (FEU)*
AMAX Auto D-Dimer	Refer to instrument application	Not specified
Trinity Biotech Auto Dimer	<150 µg/mL	Not specified
<i>ELISA quantitative assays</i>		
American Diagnostica Imuclone	0–800 ng/mL	Not specified
Zymutest D-Dimer	<400 ng/mL	>500 ng/mL
Biopool/Trinity Biotech TintElize D-Dimer	130 ng/mL	Not specified
<i>Small device or point of care</i>		
MiniQuant d-Dimer	<250 µg/mL	Same as RR *
Biosite Triage D-Dimer	0–600 ng/mL	Not specified
Vidas D-Dimer	<500 ng/mL (FEU)	>500 ng/mL (FEU)

dilution of a patient sample that still tests positive is a 1/8, then the patient's D-dimer level is greater than ($0.500 \times 8 = 4.00$ µg/mL). It is important to check product inserts for specific dilution instructions and thresholds as there are differences between the various commercially available products. Product inserts should also provide the linear limit of a product over which a "greater than" value is resulted.

Table 2
Method for qualitative slide test

Step	Details
1. Allow reagents and QC to come to room temperature.	- This usually includes antibody-coated latex particles, test cards, and mixing rods and may include buffer. Some kit reagents may not contain latex. - Allow to stand for 30 min at room temperature.
2. Label the test card for positive and negative controls and patient samples.	
3. Place the specified amount of neat patient sample on the test card in the appropriate area.	Follow manufacturer's recommendation for volume.
4. Repeat step 3 for each of the controls.	- Mix controls prior to use. - QC may need to be reconstituted prior to first use.
5. Place the appropriate amount of antibody-coated latex reagent (or other reagent) in the test area for each sample.	- The latex (or other) reagent should be adjacent to, but not mixed with, the samples. - Ensure that reagent is adequately resuspended.
6. Using separate mixing rods, combine and mix each test area individually.	
7. Rock the test card for the specified time.	- Adhere to manufacturer's mixing time and speed. - Manual or mechanical mixing may also be specified.
8. Compare the test area of the patient sample with that of the negative and positive controls.	Slide should be viewed macroscopically using a direct light source.
9. Patient samples that show agglutination contain D-dimer levels above the cutoff value.	Cutoff is specific for each commercially available kit; verify with insert and/or manufacturer.
10. Discard the test card and mixing rods in the appropriate biohazard container.	Cards are generally for single use only.

3.4 Quantitative D-Dimer Assays (Automated Analyzers)

Quantitative methods from the manufacturer for D-dimer are usually pre-set with important details (i.e., incubation times, reading wavelength, total read time etc.) in a particular instrument. For a general overview with basic instruction (see Table 3).

3.5 Interferences

As with most antibody-based lab assays, high levels of rheumatoid factor may cause interference, although some specific commercially available product inserts do state that there is no interference.

Interference by heterophile antibodies, although rare, can occur. Often the patient has an underlying malignancy such as prolymphocytic leukemia (16). This may be dealt with in a couple of

Table 3
Method for quantitative D-dimer assays (automated analyzers)

Step	Details
1. Bring reagents (usually latex suspension and buffer), controls, and calibrators to room temperature.	- Calibrators and controls and other reagents may need to be reconstituted; follow manufacturer's instruction. - Some commercial products have a calibration curve that has been predetermined for a specific lot number, and is reliable throughout the stability of that specific reagent lot.
2. Load calibrators, controls, and reagents onto the analyzer.	- Ensure that all are well mixed. - Some reagents may require mixing bars to avoid precipitation while on the analyzer (see Note 5).
3. Request a calibration to be done. A single calibration at the beginning of each lot number of reagent may be adequate, unless QC problems occur.	The calibrator is assayed using various dilutions to prepare a standard curve (use manufacturer's recommendations).
4. Request the positive and negative controls to be run.	Verify that the results of the controls are acceptable.
5. Load the patient sample and request the appropriate assays to be run (see Note 6).	- Some automated analyzers will reflex test to do a larger dilution if a patient sample is beyond the linearity of the standard curve. - This new assay dilution needs to be multiplied by the appropriate corrector factor to attain the correct result.

ways. Commercially available products to pretreat the sample removing the heterophile antibodies are used in some chemistry laboratories, and could be attempted as long as appropriate controls are assayed along with the patient sample. Also, some commercially available ELISA assays for D-dimer have a large excess of nonimmune mouse IgG in each well to minimize interference from the heterophile antibodies. Unspecified blocking agents are present in some automated methods.

Severely lipemic or hemolyzed samples may cause some interference, at least in the automated assays, as the rate of reaction is monitored as optical density or light scatter, over time. The original validation of the assay should include interference data if deemed necessary. For some assays, diluting the sample with saline, or other suitable diluent, may minimize interference (see Note 7).

The latex bead suspension for any type of assay should never be frozen. Agglutination of the suspension may occur and will be evident in the negative control.

Overmixing of the slide tests beyond the manufacturer's recommended time may result in the edges of the reaction zone drying out. The appearance of the drying slide may be judged in error to be positive.

3.6 Risk Stratification

Risk stratification, using quantitative D-dimer level as an indicator, has been investigated in a number of clinical contexts, including both acute myocardial infarction (17) and unstable angina/non-Q myocardial infarction (18). This may or may not be an important consideration, depending on the patient population being tested and the physician group interpreting the results. For dimer assays used in research studies especially, choosing an assay with a high linear limit may be useful, especially if risk stratification is considered to be a measured outcome of the study. The linear limit on some commercially available products is much higher than others. Each kit has its own defined maximum value for which a “greater than” value is assigned to all samples exceeding that value. Alternative dilutions, for some methodologies, may increase the maximum value allowed.

4 Notes

1. It is important to be careful when choosing a D-dimer assay, to consider the clinically significant cutoff value to use for that assay. This is the level above which a physician would likely expand the investigation on a patient who appeared at risk for a thrombosis. Careful review of the literature for a given product is necessary and this literature should include randomized control trials for the efficacy of the product in question. It is obviously most important to minimize false negative results, and if deemed necessary to “err on the side of caution” when choosing an appropriate cutoff value for recommendation to physicians.
2. It is very important to understand that cutoff values are strictly dependent on the assays used and are *NOT* interchangeable between various commercially available products.
3. When doing a side-by-side comparison of two or more assays, occasionally a patient will have differing results (i.e., positive or negative) for different assays, even after considering the various reference range values. The differences are likely related to the individual antibody clone(s) that is used by the manufacturer and the difference in reactivity with the various sizes of fragments. It is the author’s experience that this occurs mostly at or around the cutoff level and for patients with a higher level of dimer present, usually concordance between assays is acceptable.
4. For the ELISA-based assays, the standard precautions that apply to these assays should be followed. The plate should not dry out during the assay. Both substrate and stop solution timings should be followed accurately to ensure optimal performance.
5. Some automated latex reagents do not need to have a stirring magnet present in order to remain in suspension. These reagents also have a long onboard stability, which would make

one think a stirring magnet is essential. Latex particles are sometimes so small that they will not precipitate out from suspension, even with routine centrifugation speeds and times.

6. Automated D-dimer assays are usually very costly and can usually be resulted from a single analysis after verification of the precision of the reagent/analyzer combination in use. In actual usage, the only results that one may choose to repeat are those close to the cutoff values. Values that are normal or very high usually do not make for cost-effective repeats in most cases.
7. Lipemic samples can interfere—most frequently in automated quantitative D-dimer assays. Basically the lipemia in the sample prevents the detector from perceiving that the bead agglutination is occurring. The amount of lipemia present in a sample that would cause interference can be quantitated in order for the front-line staff to have a quick indicator when unsure if a lipemic sample has reached the threshold past which interference is a problem. By spiking samples with various amounts of lipid (much like one would do for interference testing of a chemistry assay), the point at which interference occurs can be defined. Running “lipemic indexes” on chemistry equipment or “hemoglobin blanks” on hematology equipment on these same samples converts this defined threshold into a piece of attainable data when the staff member is questioning whether or not a patient’s lipemic sample is producing interference on the particular quantitative D-dimer assay.

References

1. Hunt FA, Rylatt DB, Hart RA, Bundesen PG (1985) Serum crosslinked fibrin (XDP) and fibrinogen/fibrin degradation products (FDP) in disorders associated with activation of the coagulation of fibrinolytic systems. *BJH* 60(4):715–722
2. Adam SS, Key NS, Greenberg CS (2009) D-dimer antigen: current concepts and future prospects. *Blood* 113(13):2878–2887
3. Brown BA (1980) Coagulation. In: *Hematology*, 3rd edn. Lea & Febiger, Philadelphia, PA. p 115
4. Francis CW, Marder VJ (1994) Physiologic regulation and pathologic disorders of fibrinolysis. In: Colman RW, Hirsch J, Marder VJ (eds) *Hemostasis and thrombosis: basic principles and clinical practice*, 3rd edn. TB Lippincott Company, Philadelphia, PA, pp 1081–1083
5. Laffan M, Manning R (2006) Investigation of haemostasis. In: Lewis SM, Bain BJ, Bates I (eds) *Dacie & Lewis practical haematology*, 10th edn. Churchill Livingstone Elsevier, Edinburgh, p 386
6. Bockenstedt PL (2003) Laboratory methods in hemostasis. In: Loscalzo J, Schafer AI (eds) *Thrombosis and hemorrhage*. Lippincott Williams & Wilkins, Philadelphia, PA, pp 383–386
7. Ellman L, Carvalho A et al (1973) The Thrombo-Wellcotest as a screening test for disseminated intravascular coagulation. *New Engl J Med* 288:633
8. Dempfle CE (2006) D-Dimer assays: the current status and new assay technologies. *Thromb Res* 118:569–571
9. AMAX Accuclot D-Dimer product insert No. CRS116; revision 2003-02
10. Quantitative D-dimer for the exclusion of venous thromboembolic disease (2010) Proposed Guideline H59-P, Vol 30(No. 9)
11. Diagnostica Stago D-DI test package insert. REF 00454, 24552 02—Edition October 2003
12. Dempfle CE (2003) Evaluation of a new automated quantitative D-dimer, advanced D-dimer, in patients suspected of venous thromboembolism: a comment on correlation of D-dimer

- assays. *Blood Coagul Fibrinol* 14: 313-317
13. Walker JB, Nesheim ME (1999) The molecular weights, mass distribution, chain composition, and structure of soluble fibrin degradation products released from fibrin clot perfused with plasmin. *J Biol Chem* 274:5201-5212
 14. Pfitzner SA, Dempfle CE, Heene DL (1997) Fibrin detected in plasma of patients with disseminated intravascular coagulation by fibrin-specific antibodies consists primarily of high molecular weight factor XIIIa-crosslinked and plasmin-modified complexes partially containing fibrinopeptide A. *Thromb Haemost* 78:1067-1078
 15. Dempfle CE, Zips S, Ergul H, Heene DL, and the FACT Study Group (2001) Evaluation of 23 quantitative D-dimer assays as basis for the development of D-dimer calibrators. *Thromb Haemost* 85:671-678
 16. Hammouda W, Moroz LA (1987) False-positive reaction for fibrin degradation products due to a monoclonal (IgM lambda) cryoglobulin with warm-reactive antibody activity for rabbit IgG. *Am J Med* 82: 1263-1268
 17. Dempfle CE, Kontny F, Abildgaard U (1999) Predictive value of coagulation markers concerning clinical outcome 90 days after anterior myocardial infarction. *Thromb Haemost* 81(5):701-704
 18. Galvani M, Ferrini D, Ottani F, Eisenberg PR (1999) Early risk stratification of unstable angina/non-Q myocardial infarction: biochemical markers of coronary thrombosis. *Int J Cardiol* 68(Suppl 1):S55-S61

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