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
ATTESTATION

STATE OF FLORIDA
COUNTY OF DUVAL


On this 19th day of September 2013, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following Instruction for Use: PULSAR® II Generator Operator's Manual, document no. 70-10-1429 Rev A, presented to me by the document's custodian, Marsha Seetaram, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

United States of America
State of Florida
County of Duval

Subscribed and sworn to before me
this 19th day of September 2013

Medtronic

PULSAR® II Generator Operator's Manual

Italiano

Deutsch

Español

Nederlands

Dansk

Svenska

Suomi

Português

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Polski

Česky

Magyar

Türkçe

Norsk

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Hrvatski

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Table of Contents

Preface	iii
About This Manual	iii
Conventions Used in This Manual	iii
Help	iii
Indications	iii
Contraindications	iv
Warnings and Cautions	iv
 Introduction	 1
PULSAR II Generator	1
Electrosurgical Modes	1
Controls, Displays, and Receptacles	2
Symbols	6
Front Panel	6
Back Panel	6
 Installation	 9
Initial Inspection	9
Installation	9
Optional Wireless Footswitch	10
Preliminary Checks	10
Preliminary Functional Testing	10
 Operation	 13
Preparing for Surgery	13
Preparing for Surgery	13
Activating the Generator	15
Activation	15
Volume Adjustment	17
Preferred Set Points	18
 Cleaning and Maintenance	 19
Inspections Required Before Each Use	19
Required Annual Inspections	19
Cleaning	19
Maintenance	20
Storage	20
Environmental Protection	20

Preface

This manual and the equipment described within are for use only by qualified medical personnel possessing training in the surgical procedures to be performed. This manual is provided as a guide for using the PULSAR® II Electrosurgical Generator (Model Number: PS100-102) only.

About This Manual

The *PULSAR II Generator Operator's Manual* provides detailed information on operating and maintaining the PULSAR II Generator.

For information on accessories, refer to the appropriate instructions for use.

Conventions Used in This Manual

▲ **WARNING:** A warning indicates a hazardous condition that may result in injury or death, if not corrected or avoided.

- ◆ **CAUTION:** Alerts you to the possibility of a problem with the device associated with its use or misuse resulting in equipment damage or failure in a procedure, if not corrected or avoided.
- **IMPORTANT:** Highlights important information for a particular section.
- **NOTE:** Points out additional information that may be helpful.

Help

- 1 Read through the section of the guide specific to the operation you are performing. Refer to the table of contents and index to locate information. A glossary is included to assist you with any unfamiliar terms.
- 2 See Troubleshooting on page 21 for a list of problems and suggested solutions.
- 3 For technical support, contact 1-866-777-9400.

Indications

The PEAK Surgery System is indicated for cutting and coagulation of soft tissue during General, Plastic and Reconstructive (including but not limited to skin incisions and development of skin flaps), ENT, Gynecologic, Orthopaedic, Arthroscopic, Spinal and Neurological surgical procedures.

- Do not reuse, resterilize or reprocess “Single Use Only” labeled accessory devices, because product failure or injury may occur.
- The handpiece and Patient Return Electrodes each have unique connectors and receptacles. Do not attempt to connect to the improper receptacle because connector damage can occur, which may result in a failure. Ensure all connections are secure before using.
- Verify that the Patient Return Electrode cable is securely connected to the Patient Return Electrode receptacle. If there is an improper connection the PULSAR II Generator will not operate as indicated by the impedance sensing feature.
- Do not connect wet connectors to the Generator as this may result in equipment damage, fire at the receptacle, or injury to the patient or user.
- Do not wrap cords around metal objects as this may induce currents that could produce system performance changes, shocks, fires, or injury to the patient or surgical personnel.
- Do not use accessories other than the ones recommended in Product Accessories. Use of other accessories may result in an unintended output and/or injury to patient.
- Use active accessories with rated voltage equal to or greater than that of the Generator's maximum output voltage.

▲ Warnings for Patient Preparation

- Observe fire precautions at all times. An electrosurgical device may provide an ignition source due to sparking and heating.
- Do not use in the presence of flammable anesthetics or oxidizing gases such as nitrous oxide and oxygen. Do not activate the unit until vapors from alcohol-based skin prepping agents have dissipated. Naturally occurring gases that accumulate in body cavities can also be an ignition source.
- There is a risk of pooling of flammable solution under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before the Generator is used.
- Ensure that all oxygen circuit connections and endotracheal tubes are leak-free before and during electrosurgery use. An oxygen leak can be a possible ignition source.
- Inadvertent patient contact may result in burns. When not in use, place accessory devices in a dry and nonconductive area away from the patient.
- Position the cables for a handpiece and the Patient Return Electrode to avoid patient contact to protect against high frequency current paths to the patient, as such contact may result in patient or user injury.
- Do not allow patient contact with grounded metal objects, as such contact may result in patient or user injury.

- The active device should not be used near electrocardiograph electrodes as it can cause interference.

▲ Warnings for Active Implants

- The use of electrosurgery in the presence of internal or external active cardiac implants is potentially hazardous. Interference from the electrical current can cause active implants to malfunction. Consult the cardiac implant manufacturer for further information before proceeding with the surgery.
- To minimize the possibility of cardiac implant interference, place the Patient Return Electrode such that the electrosurgical current path is as far as possible from the pacer lead.


▲ Warnings for Use

- During use, a diminished power output may indicate failure of an electrical lead, active electrode insulation failure, or excessive eschar buildup on the active electrode tip. Do not increase the power output before checking for obvious defects or improper connections.
- Apparent low output or failure of the Generator to function correctly at the normal operating settings may indicate faulty application of the Patient Return Electrode or poor contact in its connections. In this case, the application of the neutral electrode and its connections should be checked before selecting a higher output power. Check for effective contact of the Patient Return Electrode to the patient any time that the patient is moved after initial application of the Patient Return Electrode.
- If power levels were increased to compensate diminished performance, it is recommended to reduce power to the original or a lower level upon resumption of use.
- If the system resets due to a power interruption or low voltage, the system will check for effective contact of the Patient Return Electrode, however the user should verify effective contact of the Patient Return Electrode visually prior to resuming electrosurgery.
- Do not use electrosurgery on small appendages, such as in finger surgery, as it can cause thrombosis or other unintended injury to tissue proximal to the surgical site.
- Studies have shown that smoke generated during electrosurgery may be harmful to surgical personnel. These studies recommend the use of a surgical mask and adequate ventilation of the smoke using a surgical smoke evacuator or other means.
- Observe all caution and warning notices printed on the unit.
- Operating room staff should never contact the handpiece tip while the Generator is active, as injury may result.
- The tip of a recently activated handpiece may be hot enough to cause patient burns or ignite surgical drapes or other flammable material. When not in use, store the device in an electrically insulated container or holster. Never place or rest a handpiece on the patient.
- The output power selected should be as low as possible for the intended purpose.

Introduction

PULSAR II Generator

The PULSAR II Generator is one component of the PEAK Surgery System. The Generator provides radio frequency (RF) energy to a disposable electrosurgical device. The system is used to perform electrosurgery.

PULSAR II Generator

Optional PULSAR Wireless Footswitch

Figure 1. PEAK Surgery System

Electrosurgical Modes

The PULSAR II Generator operates in the following modes:

Cut Modes

- Low Cut – Precision cutting with some hemostasis and minimal collateral damage
- Medium Cut (Pure Cut) – Precision cutting with increased hemostasis and collateral damage
- High Cut (Blend Mode) – Precision cutting with strong hemostasis


Coagulation Modes

Coagulation increases with higher settings.

- Low Coag (Pinpoint)
- High Coag (Spray)

■ IMPORTANT:

- Each cut and coagulation mode has multiple power settings available.
- Electrosurgical device models vary and may not utilize all modes. Available modes for each device are per applicable device's Instructions For Use.

Power Switch – A black rocker switch that toggles right to turn the PULSAR II Generator power on and toggles left to shut the power off.

Power Indicator – LED that illuminates green while the Generator is powered on and turns off when the Generator is powered off.

Device Connector – A seven-pin circular connector that accepts Peak PlasmaBlade devices.

Patient Return Electrode Connector – A standard two-pin connector that accepts the Patient Return Electrode connector used in procedures.

Color LCD Touch Screen – Used to prompt the user for input and to display informational and fault messages.


CUT Output Level Adjustment Buttons (on LCD) – Pressing the up button increases the Cut output level by single-integer increments. Pressing the down button decreases the Cut output level by single-integer decrements.

COAG Output Level Adjustment Buttons (on LCD) – Pressing the up button increases the Coag output level by single-integer increments. Pressing the down button decreases the Coag output level by single-integer decrements.

CUT Output Level Display – Indicates the Cut output level in single-integer increments from 1 to 10, where 1 is the lowest level and 10 is the highest level. Each device has an embedded chip that sets the Cut settings for that particular device. These settings are described in the Instructions For Use accompanying the device.

COAG Output Level Display – Indicates the Coag output level in single-integer increments. Setting 1 is the lowest level and 10 is the highest level. Each device has an embedded chip that sets the Coag settings for that particular device. These settings are described in the Instructions For Use accompanying the device.

Footswitch Connected Indicator – Indicates that the footswitch is connected to the Generator. A green footswitch icon is displayed at the top center of the display when connected.

Back Panel Layout Description

Refer to Figure 3 for a complete back panel illustration. Each control, receptacle, or panel is described in more detail below.



Figure 3. Back Panel Layout

Power Entry Module – The power entry module combines the connector for the 3-prong, Hospital-Grade power cord with removable enclosure holding two line fuses. Always use fuses of the rating shown in “Mains Input Characteristics” on page 32 of this manual.


NOTE: “Use Hospital Grade receptacle for grounding reliability” statement applies to US/North America power cords only.

Footswitch Connector – A connector that accepts a PULSAR Wireless Footswitch connector.

Equipotential Ground Connector – Standard connector for connecting common grounds.

Symbol	Description
	Equipotential ground
	Non-ionizing radiation
	Footswitch
	Explosion risk if used in the presence of flammable anesthetics
	Fuse
	Separate collection of waste at end of life as required by European Directives. Dispose of in accordance with the applicable country regulation.
	Affixed in accordance with European Council Directive 93/42/EEC.
	Authorized representative in the European community.
	In accordance with IEC 60601-1, UL 60601-1, and IEC 60601-2-2; CAN/CSA C22.2 NO. 601.1.
	Manufacturer

Installation

Initial Inspection

Unpack the Generator upon receipt and physically inspect it for any obvious damage that may have occurred during shipment. A qualified biomedical engineer or other person thoroughly familiar with electrosurgical devices should perform this inspection.

If the Generator is damaged, call 1-866-777-9400 for assistance.

Installation

The PULSAR II Generator may be placed on a stable surface.

▲ **WARNING:** The Hospital Grade power cord of the Generator should be connected to a properly polarized and grounded power source whose voltage and frequency characteristics are compatible with those listed on the nameplate of the PULSAR II Generator (located on the rear panel of the Generator). Improper grounding may result in equipment damage, fire at the receptacle, or injury to the patient or user.

NOTE: "Use Hospital Grade receptacle for grounding reliability" statement applies to US/ North America power cords only.

▲ **WARNING:** To allow for appropriate cooling, do not install the Generator in a cabinet or similar enclosure. If mounted on a shelf or near a wall, allow a 3-inch clearance around the Generator to permit free air circulation on all sides. Appropriate cooling is necessary to avoid overheating of the Generator.


Figure 5. Front Panel Showing Self-Test Screen

- The Cut activation tone and the Coag activation tone will sound to allow you to verify operation.

The screen will disappear when the self-test is complete. If the self-test is not successful, an alarm tone sounds. A Fault Code appears on the LCD screen and the Generator output is disabled. The alarm tone will sound until the Generator is turned off. In this case, call Medtronic Advanced Energy Customer Service. Do not attempt to use the Generator until the problem has been resolved.

3. If the self-test is successful, verify the following:

- If the Generator is powered on with *no electrosurgery devices connected* or with just the Patient Return Electrode connected, the Generator displays the operational screen with small dashes in the settings positions and no mode specified. The front panel showing the operational screen with no devices connected is shown in Figure 6.


Figure 6. Front Panel Showing Operational Screen with No Devices Connected

- If the Generator is powered on with an *electrosurgical device connected* (with or without the Patient Return Electrode connected), the Generator displays the screen showing the set point 1 output levels for both Cut and Coag. The front panel displaying this screen is shown in Figure 7 using a Split Foil Patient Return Electrode.

Operation

Preparing for Surgery

1. Ensure the power switch is off, and then connect the power cord to a properly grounded and polarized mating power receptacle.
2. Set the power switch to the on position. The green power ON indicator will indicate proper connection to an AC power source. The color LCD on the front panel will illuminate and show a self-test screen. As the Generator goes through its internal self diagnostics, the Generator should respond by:
 - Sounding the Cut activation tone, and then the Coag activation tone to allow verification; and
 - Replacing the self-test screen with an operational screen indicating the completion of self-test.
3. After the power-on self-tests, the Generator is ready for use.

If the Generator fails to respond as above, it has failed one of its internal tests and must not be used. An alarm tone will sound and a Fault Code will appear on the front panel LCD. The Generator output is disabled. In this case, call Medtronic Advanced Energy Customer Service at 1-866-777-9400. Do not attempt to use the Generator until the problem has been resolved.

Preparing for Surgery

1. Connect the Patient Return Electrode to the Patient Return Electrode connector on the front panel of the Generator. Ensure that the Patient Return Electrode is reliably attached to the patient's body and as close to the operating field as possible.
2. Position the Patient Return Electrode on the patient.
There are two types of Patient Return Electrodes that are compatible with the Generator: Single Foil Electrodes and Split Foil Electrodes.
 - **Single Foil Electrode** – Position the Single Foil Patient Return Electrode on the patient for the lowest impedance. Connection of the Patient Return Electrode cable to the Single Foil electrode will be indicated by a green solid bar indicator. When using Single Foil electrode types, only the cable connection is monitored by the Generator. Proper contact with the patient is *not* monitored.
 - **Split Foil Electrode** – Position the Split Foil Patient Return Electrode on the patient for the lowest impedance, as indicated by the orange bar graph. The impedance is continuously monitored and displayed on the screen. The lower the impedance, the more orange bars will be displayed. Position the Patient Return Electrode to maximize the number of orange bars and then “lock” the impedance. Impedance can be locked by either touching the bar indicator on the LCD screen or activating the Cut or Coag button on the

The user must determine the preferred settings, and set the level by pressing the UP or DOWN buttons displayed on the LCD touch screen next to the output level on the display. Output levels are displayed in single-integer increments. Disposable Device Instructions for Use will indicate the power associated with each level.

■ **IMPORTANT:** The output level settings cannot be adjusted while the Generator is activated.

4. The Generator is now ready.

Refer to “Activating the Generator” below.

Figure 8 shows an example of the front panel display when the Generator is ready for use.

▲ **WARNING:** Always keep unused, active electrodes away from the patient.


Figure 8. Front Panel when Ready for Surgery

Activating the Generator

Activation

If the Generator is not already set up, follow the set-up procedure to prepare the Generator for operation.

1. Place the electrosurgical device in the desired position in the surgical field.
2. To activate the system, depress the Cut or Coag button on the handpiece, or depress the Cut or Coag pedal on the footswitch.

Normally, the error display will disappear after 5 seconds (without re-pressing the device button), showing the normal display if the Generator detects proper Patient Return Electrode impedance. If the Generator does not detect proper impedance after the error display disappears, the bar graph will turn red as shown in Figure 11 and will remain red until proper impedance is detected.


Figure 11. Improper Locked Impedance

You may have to reapply or readjust the Patient Return Electrode in order for the Generator to detect proper Patient Return Electrode impedance. Once proper impedance is detected, the indicator will turn green.

If you are not able to adjust the Patient Return Electrode sufficiently to detect proper Patient Return Electrode impedance, unplug the patient return connector and plug it back in. This will unlock the impedance and the bar graph will turn orange. Repeat section “Preparing for Surgery” on page 13, step 2 (split foil electrode) to lock the impedance and continue with the surgery.

The Generator presents audible and visible alarms if the Generator detects improper impedance with the Patient Return Electrode and will disable Generator output. Refer to the “Patient Return Electrode Contact Quality Monitor” on page 34 for specifications related to Patient Return Electrode sensing. Refer to “Errors” on page 21 for the Error Code number displayed and other indicators while the Generator detects improper Patient Return Electrode impedance. Refer to the manufacturer’s instructions for application site and placement procedures when applying the Patient Return Electrode.

3. When the surgical procedure is completed, turn off the Generator and disconnect the electrosurgical device, footswitch, and Patient Return Electrode.

Volume Adjustment

Volume can be adjusted from the LCD touch screen’s slider next to the speaker icon.

Cleaning and Maintenance

This section contains information for ordinary cleaning and maintenance of the PULSAR II Generator. While the Generator has been designed and manufactured to high industry standards, periodic inspection and performance testing must be performed by a qualified biomedical technician for safe and effective operation.

Inspections Required Before Each Use

1. Visually inspect the Generator for physical damage. Report damage to Medtronic Advanced Energy or your biomedical department. Do not use the Generator if it is damaged.
2. Visually inspect the power cord and plug for physical damage. Replace the cord if the insulation has been breached. Do not use the Generator if the cord or plug has been damaged and has not yet been replaced.

Required Annual Inspections

1. Inspect the tightness of the power plug. If it is loose, it must be replaced with Medtronic Advanced Energy approved components.
2. Inspect the mating, cleanliness, and absence of damage to the patient connectors. Do not use the Generator if the connectors are damaged.
3. Inspect for accumulation of lint or debris within the Generator or fan vents. Do not use the Generator if lint or debris has accumulated and has not been cleared.

Cleaning

▲ WARNING: Electric shock hazard. Always unplug the Generator from the wall outlet prior to cleaning.

The Generator is not sterilizable.

1. Clean the front display, cover, and cord with a mild detergent or mild disinfecting solution and damp cloth.

It is recommended that non-flammable agents be used for cleaning and disinfection whenever possible. If flammable agents are used for cleaning, disinfecting, or as solvents, they should be allowed to evaporate before surgery.

- ◆ **CAUTION:** Do not allow fluids to enter the chassis. Do not use alcohol, caustic, corrosive, or abrasive materials on the front display, cover, and cord, as they may cause damage to the equipment. Medtronic Advanced Energy recommends following hospital procedures for cleaning the outside of the Generator after each patient.

Troubleshooting

This section identifies the possible error and fault conditions and offers common solutions for correcting malfunctions and responding to alarms.

The PULSAR II Generator has been designed and manufactured at the highest possible standards. In the rare event that the Generator fails, it has several self-diagnostic routines to aid in troubleshooting the problem. If the software detects a problem, an error code or fault code is displayed on the screen, and the Generator is disabled. The Generator will remain disabled until the detected problem is corrected. The self-diagnostic routines are only an aid for qualified technicians, and are not a substitute for evaluation of a problem by a qualified technician.

▲ WARNING: Electric shock hazard. Always unplug the Generator from the wall outlet prior to opening the cover for servicing. There are no user-serviceable parts inside the Generator. The PULSAR II Generator should only be serviced by a Medtronic Advanced Energy trained technician.

Errors

The PULSAR II Generator includes automatic self-diagnostics. If the diagnostics detect a condition that can be corrected by the operator, the system displays an error code, sounds an alarm tone, and deactivates the Generator's output. Correcting the error will enable the Generator's output.

Most error codes result from conditions in accessories attached to the Generator. The following table lists the error codes and describes the condition.

Error Code	Description
E1	Cut or Coag switch on electrosurgical device may be stuck in the ON position.
E2	Cut or Coag footpedal switch may be stuck in the ON position.
E3	Patient Return Electrode has poor connection.
E4	Generator output error.
E5	Electrosurgical device has reached end of life.
E6	Device is not recognized as a PEAK PlasmaBlade device or failed testing.
E7	Generator voltage error.
E8	Overcurrent error.

Error Code E2

Cut or Coag foot pedal switch may be stuck in the ON position.

Cause	Solution
Footswitch has Cut or Coag pedal switch stuck in the ON position when the footswitch is connected. The Error Code is displayed and the alarm will sound once.	1. Ensure that the handpiece is not in contact with the patient. 2. Disconnect and reconnect the footswitch. 3. If the Error Code reappears, turn the Generator power OFF and then turn the power ON again. 4. If the problem persists, replace the footswitch and repeat this procedure. 5. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Error Code E3

Patient Return Electrode has poor connection.

Cause	Solution
The Error Code is displayed when the impedance is out of range. The Error Code is displayed for 5 seconds and the alarm will sound once. The Patient Return Electrode impedance icon will turn red and stay red as long as the impedance is out of range. The Patient Return Electrode impedance is continuously checked.	1. Follow the manufacturer's instructions for Patient Return Electrode placement. 2. Verify that the Patient Return Electrode is connected to the Generator and is making proper contact with the patient. 3. If the impedance icon remains red, replace the Patient Return Electrode and connector. 4. Repeat preparations for surgery. Position the Patient Return Electrode making contact for lowest impedance. 5. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Error Code E7

Generator voltage error.

Cause	Solution
Measurements of the output voltage are being made while the Generator has active output. If these measurements are outside the normal range, the output will be disabled. The Error Code is displayed and the alarm will sound once.	1. Pressing the footswitch or handpiece button again will clear this screen and reactivate output if the error is no longer detected. The handpiece should not be in contact with the patient while this is done. 2. If the Error Code reappears, replace the handpiece. 3. If the Error Code reappears, disconnect the handpiece and the Patient Return Electrode. Turn the Generator power OFF and then turn the power ON again. Reconnect the handpiece and Patient Return Electrode and try again. 4. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Error Code E8

Overcurrent error.

Cause	Solution
Measurements of the output current are being made while the Generator has active output. If these measurements are over the normal range, the output will be disabled. The Error Code is displayed and the alarm will sound once.	1. Pressing the footswitch or handpiece button again will clear this screen and reactivate output if the error is no longer detected. The handpiece should not be in contact with the patient while this is done. 2. If the Error Code reappears, replace the handpiece. 3. If the Error Code reappears, disconnect the handpiece and the Patient Return Electrode. Turn the Generator power OFF and then turn the power ON again. Reconnect the handpiece and Patient Return Electrode and try again. 4. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Fault Code F1

Internal temperature of the Generator has exceeded the limit.

Cause	Solution
Internal temperature exceeded limit while turned ON. The Fault Code is displayed and the alarm will sound once.	1. Leave the Generator ON with the internal fan running to cool the Generator. Do not use the Generator until it has cooled sufficiently. Once the Generator has cooled, the Fault Code will no longer be displayed. 2. If the Fault Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Fault Code F2

Self-test fault.

Cause	Solution
Self-test diagnostics are run when the Generator is powered ON. This Fault Code indicates that one or more of the self-tests have failed. The individual self-test that failed will be recorded in the Generator's data flash for retrieval during maintenance. The Fault Code is displayed and the alarm will sound continuously.	1. Turn the Generator power OFF and then turn the power ON again. Do not press buttons or activate handpieces during the self-test. 2. If the Fault Code reappears, disconnect the handpiece and Patient Return Electrode. Turn the Generator power OFF and then turn the power ON again. 3. If the Fault Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Fault Code F3

RF Generator not calibrated or calibration was lost.

Cause	Solution
The RF Module has detected that it is not calibrated or calibration values have been lost. The Fault Code is displayed and the alarm will sound continuously until the Generator is turned off.	1. Turn the Generator power OFF and then turn the power ON again. Reconnect the handpiece and Patient Return Electrode and try again. 2. If the Fault Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Glossary

This section provides definitions of terms relating to electrosurgery. The information below is provided as a convenience and should not be interpreted as an authoritative discussion of the subject.

active electrode	The electrode or electrode assembly at which the electrosurgical effect is intended. It usually has a small contact surface area and provides a high current density to achieve the desired surgical effect.
blend	A modified cut waveform incorporating some off-time, thus allowing tissue to cool, and providing varying degrees of hemostasis.
coagulation	The clotting of blood or dehydration of tissue with no cutting effect. Electrosurgical coagulation often incorporates intermittent bursts of high-voltage, low-current electricity.
cutting	Occurs when intense sparks focus the energy and vaporize the tissue. It results from high current density in the tissue causing cellular fluid to burst into steam and disrupt the structure. Voltage is low and current flow is high. In order to induce a cutting effect in human tissue, voltage must exceed 200 Vp.
desiccation	A procedure where a small amount of surface tissue is dried out by placing the active electrode in contact with the tissue. Desiccation differs from fulguration in that peak voltage is lower resulting in the inability of the current to arc through air to tissue. Direct contact between electrode and tissue is required. See fulguration.
dispersive electrode	The electrode at which no electrosurgical effect is intended or desired. It is usually large in area in order to provide a low current density so that no electrosurgical effect occurs at that site. It is also known as a patient plate, patient pad, plate electrode, return electrode, neutral electrode, inactive electrode. It is sometimes (inaccurately) referred to as a grounding plate.
disposable accessory	An electrosurgical accessory, such as active electrodes, handles, dispersive electrodes, etc. It is not intended to be used more than once.

Specifications

PULSAR II Generator

All specifications are nominal and are within $\pm 20\%$ of a stated value at room temperature (77°F/ 25°C) and a nominal input power voltage. All specifications are subject to change without notice.

General

Classification:	IEC 60601-1, Class I
Type:	CF
Internal Design:	Solid state
Output Configuration:	Isolated (RF Floating)
Cooling:	Forced air (fan)
Designed to Meet:	UL60601-1; IEC 60601-1, 60601-1-4, 60601-1-2, 60601-2-2; CAN/CSA C22.2 NO. 601.1.

Dimensions and Weight

Width:	13.8 in (35 cm)
Length:	15.2 in (38.5 cm)
Height:	7.1 in (18 cm)
Weight:	<26.5 lb (12 kg)

Operating Parameters

Temperature:	10°C to 40°C
Humidity:	30% to 75% relative humidity, non-condensing
Degree of Protection Against Ingress of Water:	IPX0

Storage Parameters

Temperature:	-20°C to 60°C
Humidity:	Non-condensing relative humidity
Shipping Conditions for Packaged Unit:	-40°C to 70°C; 10–100% relative humidity, including condensation

Maximum RF Voltage Output

Mode	Maximum Voltage
Low Cut 1	1600
Low Cut 2	1200
Medium Cut	1600
High Cut (Blend 1)	2200
High Cut (Blend 2)	3100
High Cut (Blend 3)	3900
Low Coag 1	3500
Low Coag 2	4600
High Coag	5500

High Frequency (RF) Leakage Current

- RF leakage current <150 mA

Output Characteristics

Mode	Power (watts) ±20%	Rated Load (ohms)	Operating Frequency (Typical, kHz)	Duty Cycle (Typical, %)
Low Cut 1	40	100	1000	3.75
Low Cut 2	160	200	1000	100
Medium Cut	60	400	1000	100
High Cut (Blend 1)	60	400	1000	75
High Cut (Blend 2)	60	400	1000	50
High Cut (Blend 3)	80	400	1000	100
Low Coag 1	60	400	1000	19
Low Coag 2	80	400	1000	100
High Coag	60	500	1000	12

User Duty Cycle

- Cut Mode: 10 seconds ON, 30 seconds OFF
- Coag Mode: 10 seconds ON, 30 seconds OFF

Guidance and Manufacturer's Declaration - Electromagnetic Immunity		
The PULSAR II Generator is intended for use in the electromagnetic environment specified below. The customer or user of the PULSAR II Generator should assure that it is used in such an environment.		
Immunity Test	Test Level and Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) EN 61000-4-2	$\pm 2, 4, 6$ kV Contact discharge $\pm 2, 4, 8$ kV Air discharge	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated Immunity EN 61000-4-3	80 Mhz to 2.5 GHz 3 V/m 80% @ 1 kHz	
Conducted Immunity (AC Power) (I/O Lines) EN 61000-4-6	0.15 to 80 MHz 3 Vrms 1 kHz AC Mains	
Electrical fast transients (AC Power) (I/O Lines) EN 61000-4-4	± 2 kV AC Mains ± 1 kV for I/O Lines 5/50 5 kHz	Mains power quality should be that of a typical commercial or hospital environment.
Surge Line to Line (AC Power) IEC 61000-4-5	± 1 kV Line to line ± 2 kV Line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Magnetic Immunity EN 61000-4-8	3 A/m 50/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips and interruptions EN 61000-4-11	>95% dip in U_t for 0.5 cycle 60% dip in U_t for 5 cycles 30% dip in U_t for 25 cycles >95% dip in U_t for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the PULSAR II Generator requires continued operation during power mains interruptions, it is recommended that the PULSAR II Generator be powered from an uninterruptible power supply or a battery.
NOTE: U_T is the AC mains voltage prior to application of the test level.		

Other Specifications

- Power Cord:** Generators have a power entry connector to allow connection with a US Hospital Grade power cord or an EU plug power cord. Units accept 100–240 VAC at a frequency of 50/60 Hz, fused at 6.3 amps.
- Display Size:** 6.4 in diagonal (3.8 in tall x 5.2 in wide)


Figure 16. Output Power vs Impedance for Medium Cut


Figure 17. Output Power vs Impedance for High Cut (Blend 1)


Figure 20. Output Power vs Impedance for Low Coag 1


Figure 21. Output Power vs Impedance for Low Coag 2

Limited Warranty

Medtronic Advanced Energy LLC warrants the PULSAR II Generator against defects in materials or workmanship under normal use and preventive maintenance for the respective warranty periods shown below. Medtronic Advanced Energy will repair or replace the unit without charge, at its option, of any product, or part thereof, which has been returned to Medtronic Advanced Energy LLC or its Distributor within the applicable time period shown below. This warranty does not apply to any product or part that has been adversely affected due to use with devices manufactured or distributed by parties not authorized by Medtronic Advanced Energy LLC; repaired or altered outside Medtronic Advanced Energy's factory; subjected to improper use, negligence or accident or not used in accordance with the instructions provided in this manual. Standard hospital preventive maintenance procedures should be followed and performed by qualified biomedical service personnel and are not covered in the warranty.

- Medtronic Advanced Energy's products are warranted for the following periods after delivery to the original purchaser:
- PULSAR II Generator - One (1) Year, Parts and Labor
- PULSAR Wireless Footswitch - One (1) Year, Parts and Labor

Product Accessories

The following accessories are compatible for use with the PULSAR II Generator. Only Medtronic Advanced Energy recommended accessories should be used with the Generator.

- PEAK PlasmaBlade Devices
- PULSAR Wireless Footswitch
- Standard Patient Return Electrode or Pad
 - Single Foil Electrodes – Valleylab Standard PolyHesive™ II E7506
 - Split Foil Electrodes – Adult REM Valleylab PolyHesive II E7507

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Output Characteristic Curves


Figure 14. Output Power vs Impedance for Low Cut 1


Figure 15. Output Power vs Impedance for Low Cut 2

Patient Return Electrode Contact Quality Monitor

The system presents audible and visible alarms when it senses no Patient Return Electrode.

Cord Fault Single Foil: Continuous measurement: 20 ± 5 Ohms in resistance will cause an alarm. When the alarm condition exists, the system deactivates output power.

Split Foil: Continuous measurement: Once the system establishes the split foil electrode impedance, an increase of 30% or more will cause an alarm. When the alarm condition exists, the system deactivates output power.

Electromagnetic Compatibility

The PULSAR II Generator conforms to electromagnetic compatibility standard EN/IEC 60601-1-2.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The PULSAR II Generator is intended for use in the electromagnetic environment specified below. The customer or user of the PULSAR II Generator should assure that it is used in such an environment.		
Test Type	Test Level/ Compliance Level	Electromagnetic Environment – Guidance
Conducted emissions EN55011:2007; FCC Part 18:2009	Class A Group 1 30 MHz to 1 GHz	The PULSAR II Generator must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Radiated emissions EN55011:2007; FCC Part 18:2009	Class A Group 1 30 MHz to 1 GHz	The PULSAR II Generator is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded. ▲ WARNING: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the PULSAR II Generator or shielding the location.
Harmonics EN 61000-3-2	Class A Device	
Flicker EN 61000-3-3	Class A Device	

Audio Specifications

The audio levels stated below are for activation tones (Cut and Coag) and alarm tones (return electrode and system alarms) at a distance of one meter. Alarm tones meet the requirements for IEC 60601-2-2.

	Activation Tone	Alarm Tone
Volume:	45 to 65 dB, adjustable	70 $\pm$ 5 dB, not adjustable
Frequency:	Cut: 1.5 kHz Coag: 600 Hz	2 kHz for 0.5 seconds 1 kHz for 0.5 seconds
Duration:	Continuous while the Generator is activated	2 seconds for Errors, continuous for Faults

Mains Input Characteristics

Frequency Hz	Mains Voltage (V)			Mains Current (A)	Fuses*
	Nominal	Minimum	Maximum	Maximum	
50/60	115/230	100	240	5.3	T6.3A

* Fuses are Type T, High Breaking

Low Frequency (50/60Hz) Leakage Current

- Earth leakage current in normal condition <300 μ A
- Earth leakage current in single-fault condition <1 mA
- Enclosure (touch) leakage current in normal condition <100 μ A
- Enclosure (touch) leakage current in single fault condition <500 μ A
- Patient leakage current in normal condition <10 μ A
- Patient leakage current in single fault condition <50 μ A

Maximum RF Power Output

Mode	Maximum Power
Low Cut	160 watts into 200 ohms
Medium Cut	60 watts into 400 ohms
High Cut	80 watts into 400 ohms
Low Coag	80 watts into 400 ohms
High Coag	60 watts into 500 ohms

duty cycle	The proportion of time (expressed in percentage) that a current or device is on versus off. Duty cycle may be used when referring to current waveforms that are repetitive. Thus high-frequency current would be on using a shorter period of time than a low-frequency current. Duty cycle is also used in reference to electrical components or equipment. For example, some equipment is designed to be used continuously, that is, with a duty cycle rating of 100%, while other equipment may be rated for intermittent use, that is less than 100% duty cycle. Most ESU's are designed to be used intermittently. Typically, they are rated for 50% duty cycles. Use of any equipment beyond it's duty cycle rating may result in premature failure.
electrosurgery (surgical diathermy)	The generation and delivery of radio frequency current between an active electrode and a dispersive electrode for the purposes of dehydration of tissue. Electrosurgery also includes cutting or vaporizing (tissue explosion). In contrast to electrocautery, the electric current actually passes through the patient.
electrosurgical accessory	Equipment used in conjunction with the electrosurgical generator to accomplish electrosurgery. These include, but are not limited to footswitch, cable, dispersive electrode, and active electrode.
fulguration (spray technique)	Coagulating tissue or blood by means of radio frequency electric arcs. In contrast to desiccation, the active electrode is not in good electrical contact with the tissue, and arcs jump from the electrode to the tissue.
impedance (measured in ohms)	Total opposition, both resistive and reactive, a circuit offers to the flow of alternating current at a given frequency.
mode (operating)	Each of the distinct ways in which the electrosurgical unit can be operated with electrosurgical output, for example, cutting, coagulation, blended, spray (fulguration).
monopolar	The traditional form of electrosurgical circuit, which uses an active electrode to apply the therapeutic current to the surgical site, and a dispersive electrode (patient plate) to return the current to the ESU.
non-volatile memory	Memory type that can retain the stored information even when not powered.
power	The rate at which energy is produced or consumed. Power is equal to voltage times current, or resistance times current squared. The unit of measure of power is the watt.

Non-Recoverable Faults

The PULSAR II Generator includes automatic self-diagnostics. If the diagnostics detect a non-recoverable fault, the system displays a fault code, sounds an audible tone, and deactivates the Generator's output. All non-recoverable faults indicate a possible problem with the equipment. Power must be turned off and then back on before operation can continue. If the Self-Test does not successfully complete, refer the equipment to Service Personnel.

Most non-recoverable fault codes result from faults in the Generator. The following table lists the non-recoverable fault codes.

Non-Recoverable Fault Code	Description
F1	Internal temperature of the Generator has exceeded the limit.
F2	Self-test fault.
F3	RF Generator not calibrated or calibration was lost.

Fault Code Details

A detailed description of the cause of each fault code is provided below, along with steps to follow. Figure 13 shows an example of a Fault screen.


Figure 13. Fault Code F1

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Error Code E5

Electrosurgical device has reached end of life.

Cause	Solution
For safety reasons, each device has a limited active life span. If the time is exceeded, then the Error Code is displayed and the alarm will sound once.	1. Replace the electrosurgical device. 2. The Error Code screen will disappear if the device is recognized, tested, and has a life expectancy. 3. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Error Code E6

Device is not recognized as a PEAK PlasmaBlade device or has failed testing.

Cause	Solution
The device is connected to the Generator and it is not recognized as a PEAK PlasmaBlade device or it has failed testing. The Error Code is displayed and the alarm will sound once.	1. Replace the device. The Error Code screen will disappear if the device is recognized, tested, and has a life expectancy. 2. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Error Code Details

A detailed description of the cause of each error is provided below, along with steps to follow. Figure 12 shows an example of an Error screen.


Figure 12. Error Code E1

Error Code E1

Cut or Coag switch on electrosurgical device may be stuck in ON position.

Cause	Solution
Electrosurgical device has Cut or Coag button stuck in the ON position when the device is connected. The Error Code is displayed and the alarm will sound once.	1. Ensure that the handpiece is not in contact with the patient. 2. Disconnect and reconnect the handpiece. 3. If the Error Code reappears, turn the Generator power OFF and then turn the power ON again. 4. If the problem persists, replace the handpiece and repeat this procedure. 5. If the Error Code reappears, record the number and contact Medtronic Advanced Energy Customer Service. Do not use the Generator.

Refer to the *PULSAR Wireless Footswitch Operator's Manual* for detailed cleaning instructions for the footswitch and receiver.

Maintenance

The PULSAR II Generator must be performance tested by a hospital qualified biomedical technician at least every year. Follow your hospital's procedures for periodic performance verification of electrosurgical generators. Medtronic Advanced Energy recommends the following checks:

- Power Output verification for each mode of operation
- RF Leakage
- AC Leakage
- AC Ground Integrity

Storage

The Generator must be thoroughly checked by a qualified biomedical engineer if it has been stored for longer than 6 months.

Allow the Generator to remain at room temperature for at least 1 hour if it has been stored at extreme temperatures. Refer to "Storage Parameters" on page 31 for the limits for storing the PULSAR II Generator.

Environmental Protection

Retain the shipping container and packing material in the event that the Generator needs to be returned for repair or service. At the end of the equipment's life, dispose of it in accordance with your local regulations. The materials in the Generator include aluminum, steel, copper, thermoplastics, and electronic components.

Preferred Set Points

There are four memory set points. Set points 1, 2, and 3 can be programmed to save the preferred CUT and COAG settings and the volume level. Memory set point 4 is reserved for the default CUT and COAG settings and volume level.

- To save the current CUT and COAG settings and volume level, press and hold the memory set point for 2 seconds.
- To retrieve the preferred memory settings, press the memory set point.

While the Generator is activated, the appropriate audible tone (Cut or Coag) is sounded, the yellow background will brighten during Cut output, and the blue background will brighten during Coag output, as shown in Figure 9 and Figure 10. Releasing the button or pedal will deactivate the Generator and the display will show a dimly lit background.


Figure 9. Cut Activation


Figure 10. Coag Activation

The Patient Return Electrode green solid bar indicator shows connection to a Single Foil electrode. When the solid bar is green, the Generator confirms that the cable is plugged in and an electrode is attached. The user must confirm proper electrode attachment to the patient. When the Generator detects increased impedance, the solid bar turns red, an alarm will sound, and output will be disabled. If the Patient Return Electrode becomes disconnected, the Generator will alarm and the indicator will disappear.

The Patient Return Electrode green bar graph will show fewer and fewer bars as the impedance increases, indicating poorer Patient Return Electrode connection when using a Split Foil electrode. When the Generator detects that the impedance has increased by 30% or more, an audible and visible error alarm will occur and output will be disabled. Pressing the Electrosurgical device button while the error is displayed will show the normal display and output will be activated if the Generator detects proper Patient Return Electrode impedance. Re-pressing the device button while the Generator detects improper impedance causes the error alarm to reoccur.

electrosurgical device or footswitch. Once locked, the bar graph color will change from orange to green and all bars should be displayed.

■ **IMPORTANT:** This feature is only available for Split Foil Patient Return Electrodes.

The Generator presents audible and visible alarms if it detects improper impedance with the Patient Return Electrode and will disable Generator output. Refer to the "Patient Return Electrode Contact Quality Monitor" on page 34 for specifications related to Patient Return Electrode sensing. Refer to the "Errors" on page 21 for the Error Code number displayed and other indicators when the Generator detects improper Patient Return Electrode impedance. Alarms will be presented only for the compatible Patient Return Electrodes specified in "Product Accessories" on page 43.

▲ **WARNING:** Refer to the manufacturer's instructions for application site and placement procedures when applying the Patient Return Electrode.

▲ **WARNING:** Do not depend solely on the impedance indicator for confirmation of good Patient Return Electrode application. Qualified personnel should make the final decision on proper Patient Return Electrode placement.

3. Plug the PEAK PlasmaBlade device into the electrosurgical device connector.

Upon connection, the Generator will alarm if the device is not recognized as a PEAK PlasmaBlade device, if it fails a unit test, or if it has no life expectancy. The display will appear as shown in Figure 7. It will also display the default output levels for Cut and Coag. The green indicator in the upper-right corner of the display shows the Patient Return Electrode impedance when using a Split Foil Patient Return Electrode.

▲ **WARNING:** Always stow unused devices in an electrically insulated location such as a safety holster.

The optional footswitch may be used. If the footswitch is used, plug the wireless receiver into the footswitch connector. While using a footswitch, the output will be delivered via the electrosurgical device. The activation button on the device will continue to function while a footswitch is connected to the Generator. When the footswitch is activated, the display background will change from:

- white to yellow when the Cut pedal is pushed.
- white to blue when the Coag pedal is pushed.

■ **IMPORTANT:** Use only the PULSAR Wireless Footswitch.

■ **IMPORTANT:** Electrosurgical devices require a Patient Return Electrode. The Generator must detect proper Patient Return Electrode impedance before the Generator output can be active.


Figure 7. Front Panel Displaying the Screen with Electrosurgical Device Connected

Optional Wireless Footswitch

1. If using the optional footswitch, plug the wireless receiver into the footswitch connector of the PULSAR II Generator (see Figure 3 on page 5).
2. Turn on the PULSAR II Generator and ensure that the wireless receiver's power-on LED is green, indicating that the receiver is powered on.


Figure 4. Wireless Receiver

▲ WARNING: Only the primary surgeon using the handpiece should operate the footswitch. Unintentional activation may occur if the footswitch is activated by a separate user, which may result in patient or user injury.

Refer to the *PULSAR Wireless Footswitch Operator's Manual* for detailed product information.

Preliminary Checks

Prior to installing and using the Generator, it is strongly recommended that a qualified maintenance technician ensure proper and safe operation by testing the performance of the Generator. If problems are observed, the Generator should not be used until the problems have been resolved by a qualified technician. Call Medtronic Advanced Energy Customer Service at 1-866-777-9400.

Preliminary Functional Testing

Preliminary functional testing should be conducted by a qualified technician or operating room personnel. Personnel conducting this testing should have fully read and understood this manual.

1. Turn on the Generator by pressing the power switch.
2. Verify the following:
 - A self-test screen appears indicating that the system is performing a power-on self-test as shown in Figure 5. The software version is also shown at the bottom of the screen. If the self-test screen does not appear, cycle the power and try again. If the same problem continues, call Medtronic Advanced Energy Customer Service at 1-866-777-9400. Do not attempt to use the Generator until the problem has been resolved.

Symbols

The following symbols are found on the front and back panels of the Generator.

Front Panel

Symbol	Description
	Defibrillation-proof type CF equipment
	RF Isolated - Patient connections are isolated from earth at high frequency
	Warning: Dangerous voltage
	Caution, consult accompanying documents
	Power ON - To indicate connection to the mains
	Power OFF - To indicate disconnection from the mains
	Consult instructions for use

Back Panel

Symbol	Description
	Model number
	Serial number
	Date of manufacture
	Caution, consult accompanying documents

Patient Return Electrode Indicators:

Single Foil Electrode Indicator – Indicates connection to a Single Foil Patient Return Electrode. The solid bar will turn green when a Single Foil electrode is properly connected to the Generator. The green bar will turn red when the Generator detects an improper cable connection to the electrode. It does *not* monitor contact with the patient. Use a Split Foil electrode for contact impedance monitoring. If the bar turns red, the Patient Return Electrode should be inspected and repositioned or replaced to ensure appropriate placement and connection.

Split Foil Electrode Indicator – Indicates the impedance level of the Patient Return Electrode. The bar graph will remain green with proper Patient Return Electrode impedance. The bar graph will turn red when the Generator does not detect proper Patient Return Electrode impedance. If the bar graph turns red, the Patient Return Electrode should be inspected and repositioned or replaced to ensure appropriate placement and connection.

Preset Memory Buttons – Pressing and holding the preset memory button for 2 seconds will save the current Cut and Coag set points, as well as the volume level, into the non-volatile memory. Briefly pressing the preset memory button will recall saved Cut and Coag set points and the volume level and display them on the screen.

Volume Control Slider – Moving the volume control slider nub to the left will decrease the volume, while moving it to the right will increase the volume.

Other features of the PULSAR II Generator include:

- Proprietary devices developed by Medtronic Advanced Energy to deliver RF energy to the patient. Only PEAK PlasmaBlade devices may be used with the PULSAR II Generator.
- Device hand control
- Device footswitch control (infrared wireless footswitch)
- Four preset memory buttons to retain and recall user preferences (default set points and volume level)
- Connection to Patient Return Electrode and impedance sensing feature
- Color LCD with touch screen feature to prompt user for input and to display informational, error, and fault messages
- Audio feedback—activation and alarm tones

Controls, Displays, and Receptacles

This section describes each component of the PULSAR II Generator and its function. The controls, displays, and receptacle for electrosurgical devices, Patient Return Electrodes, and the footswitch are located on the front and rear panels of the Generator.

▲ WARNING: Read all warnings and instructions provided with the Generator prior to use.

Front Panel Layout Description

Refer to Figure 2 for a complete front panel illustration. Each display, control, or receptacle is described in more detail below.


Figure 2. Front Panel Layout

- Failure of the high-frequency surgical equipment could result in an unintended increase of output power.
- If neuromuscular stimulation is observed, check all connections to the Generator, Patient Return Electrode, and active electrodes to look for a possible metal-to-metal spark. Stimulation can be alleviated by using low-power settings or low-voltage coag. If the problem continues, the Generator should be checked by qualified service personnel for abnormal current leakage.

▲ Warnings for Testing or Servicing

- Do not remove the Generator cover due to electrical shock hazard. There are no user serviceable parts inside.
- Do not modify any component of the Generator.
- Avoid contact with the output leads when the Generator is activated as this may result in injury. Periodically inspect the test leads used for the output connections for obvious defects.
- The Generator is not designed to operate for extended periods of continuous output. When testing, it is recommended that duty cycles be limited to 25% with maximum activation times of 10 seconds into a load greater than or equal to 600 ohms. Use for an extended period of time may result in overheating and equipment damage.

◆ Cautions

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- The PULSAR II Generator should only be used by qualified medical personnel possessing training in the surgical procedures to be performed.
- Prior to initial use of the PEAK Surgery System, ensure that all package inserts, including warnings, cautions, and instructions for use are read and understood.
- Use the lowest power setting possible and the shortest activation time possible to achieve the desired end effect.
- Periodic maintenance should be performed by a hospital qualified biomedical technician or by a qualified Medtronic representative.
- Minimizing operating temperature and extreme thermal cycles will extend the equipment life.
- The heat dissipation capability of the Generator heat sink may be severely impaired by activating the PULSAR II Generator in other than a normal operating position. Testing or using the unit in any other position should be avoided.
- The PULSAR II Generator is intended for use only with a PEAK PlasmaBlade single-use disposable device.

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

ATTESTATION

STATE OF FLORIDA
COUNTY OF DUVAL


On this 19th day of September 2013, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following Instruction for Use: PULSAR® II Wireless Footswitch Operators Manual, document no. LBL-00098 Rev C, presented to me by the document's custodian, Marsha Seetaram, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
Medtronic

PULSAR® Wireless Footswitch Operators Manual

P/N LBL-00098 Rev C
Page 1 of 257

Contents

English

Safety	9
Definitions of Terminology	10
Electromagnetic Interference	12
Equipment Symbols	13
Channel/Coding	14
Receiver Front/Rear View	15
Transmitter	17
Battery Replacement	18
Cleaning	19
Specifications	20

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Safety

The information presented in this section is important for the safety of both the patient and the operator and also serves to enhance equipment reliability. This section describes how the terms Warning, Caution, Important, and Note are used throughout the manual.

General Information

General Use

If the footswitch is cold to the touch or below ambient temperature, allow it to stabilize before use.

The use of ACCESSORY equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of safety of the resulting system. Consideration relating to the choice shall include:

- use of the accessory in the PATIENT VICINITY
- evidence that the safety certification of the ACCESSORY has been performed in accordance to the appropriate IEC 60601-1 and/or IEC 60601-1-1 harmonized national standard.

Periodically, whenever the integrity of the switch is in doubt, test all functions.

Warnings:

- Explosion hazard-** Do not use this equipment in the presence of flammable anesthetics.
- Interfacing to-Equipment** Footswitches must be interfaced with other equipment. Be certain to consult manufacturer's specifications to maintain safe operation.

Cautions:

- Annual servicing-** For continued safety and performance of the switch, it is recommended that the functionality and electrical safety of the switch be verified on an annual basis.
- Daily testing-** It is essential that the footswitch be inspected everyday or before use.
- Performance-** The infrared wireless footswitch system operates on a wavelength of 870nm. It should be tested to assure compatibility with any device it is connected to or environment it is working in. Report all problems experienced with the footswitch. If the footswitch is not working properly, contact your service representative for service. The footswitch should not be used if it is not working properly.
- Rx Only-** In the United States of America, Federal Law restricts this device to sale by or on the order of a Physician.


Equipment Symbols

The following is a list of symbols used on the product.

Table 1-2. Equipment Symbols


	ATTENTION: Consult accompanying documents
	POWER ON
	LOW BATTERY
	Follow instruction for use.
	Date of manufacture
	Separate collection of waste at end of life as required by European Directives. Dispose of in accordance with the applicable country regulation.
	Manufacturer

Receiver Front View

	Name	Description
A	Power On indicator.	The power On indicator lights green to indicate the unit is powered
B	Switch function Indicators 1-5	LED's will light when switch(s) are actuated.
C	Low battery indicator foot switch	Red LED will flash when the batteries are close to the end of their life. The LED will flash for 3-4 hours. When LED stays on constantly there is 3-4 hours of battery life remaining. However batteries should be replaced immediately.
D	IR Key label	This label indicates the manufacturer's code and channel that the receiver uses to pair with a footswitch. This label is color coded and the footswitch has an identical label (page 17). Refer to page 14 for details of channel coding.

Transmitter

	Description
A	Primary switches Single or Dual Function
B	IR Key label. This label indicates the manufacturer's code and channel that the receiver uses to pair with a footswitch. This label is color coded and the receiver has an identical label (page 15). Refer to page 14 for details of channel coding.

Cleaning

The following cleaning instructions are provided.

CAUTIONS:

ABRASION - Do not use abrasive cloth, sharp objects, or abrasive cleaners.

DISCONNECTION - Detach the interconnect cables from the receiver.

IMMERSION - Do not immerse the receiver, cables or connectors under running water.

Transmitter

The transmitter is *IP68 rated* and *can be* completely immersed briefly. Dampen a non-abrasive cloth with one of the following products; then wring out until slightly wet and gently rub soiled area until clean.

- Isopropyl alcohol
- Soap and water
- Cidex
- Sodium Hypochlorite 5.25% (Bleach) diluted 10:1

Receiver

CAUTION: The receiver is *IPX1* rated and *cannot* be immersed.

Dampen a non-abrasive cloth with one of the following products; then wring out until slightly wet and gently rub soiled area until clean.

- Isopropyl alcohol
- Soap and water
- Cidex
- Bleach Sodium Hypochlorite 5.25% (Bleach) diluted 10:1

Wipe any fluids from the surface of the receiver.

This page intentionally left blank.

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Copyright