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Important information

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favour of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper use, cleaning, disinfection and, if necessary, sterilization of the intubation videoscope. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

WARNING: KARL STORZ instruments are not sterile when delivered, and must be cleaned, disinfected and, if necessary, sterilized prior to initial use and each subsequent use.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures, and descriptions in this manual.

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
7	Leakage test28
7.1	Performing the leakage test – dry28
7.2	Performing the leakage test – in liquid29
8	Disassembly.....30
9	Manual wipe-down disinfection.....30
9.1	Manual reprocessing.....30
9.1.1	Manual cleaning.....30
9.1.2	Manual disinfection31
9.1.3	Manual High Level Disinfection (HLD).....31
9.2	Drying31
10	Machine reprocessing.....32
10.1	Connecting32
10.2	Machine cleaning/chemical-thermal disinfection32
11	Assembly, inspection and care.....33
12	Packaging.....34
13	Sterilization34
13.1	From SN PR 1191 onwards: Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®35
13.2	From SN PR 1191 onwards: Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® V-PRO™35
13.3	Ethylene oxide sterilization (EO).....35
14	Limits of reprocessing36
15	Storage36
15.1	Storage.....37
16	Reprocessing the tube holder.....38
16.1	General warnings.....38
17	Precleaning directly after use.....39
18	Disassembly.....39
19	Manual reprocessing.....39
19.1	Manual cleaning.....39
19.1.1	Manual disinfection40
19.1.2	Manual High Level Disinfection (HLD).....40
19.2	Drying40

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20	Machine reprocessing	41
20.1	Connecting	41
20.1.1	Machine cleaning/ chemical-thermal disinfection.....	41
20.1.2	Machine cleaning/thermal disinfection	41
21	Assembly, inspection and care	42
22	Packaging	42
23	Sterilization	42
23.1	Steam sterilization using the fractionated prevacuum procedure.....	43
23.2	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	43
23.3	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® V-PRO™	44
23.4	Ethylene oxide sterilization (EO).....	44
23.5	Limits of reprocessing.....	44
24	Technical description	45
24.1	System description	45
24.2	Technical data.....	46
24.2.1	C-MAC VS 10331 BX	46
25	Service, repair and disposal	47
25.1	Servicing and repair	47
25.2	Maintenance	47
25.3	Limitation of liability	47
25.4	Warranty	47
25.5	Electromagnetic Compatibility.....	48
25.6	Directive compliance.....	48
25.7	Important information on dispatch.....	48
25.8	Disposal.....	49
25.9	Troubleshooting	50
26	Appendix	52
26.1	Accessories	52
27	Subsidiaries	56

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


**Controls, displays,
connectors and their uses**

- ① Connection for leakage tester or pressure compensation cap (vent port)
- ② Handle
- ③ Socket for video connecting cable ②⑥
- ④ Multifunction button (BlueButton)
- ⑤ Flexible sheath for deflection
- ⑥ Distal tip
- ⑦ Deflection lever
- ⑧ Adaptor for tube fixation
- ⑨ O₂ insufflation connection
- ⑩ Leakage tester 13242 XL
- ⑪ Pressure compensation cap 11025 E
- ⑫ Monitor
- ⑬ Secondary connection for videoscope
- ⑭ Primary connection for videoscope
- ⑮ Mains connection
- ⑯ Reset button
- ⑰ Holder for SD memory card
- ⑱ Status bar
- ⑲ Pilot lamp
- ㉑ On/Off switch
- ㉑ Charging indicator lamp (orange)
- ㉒ Soft keys
- ㉓ Function symbols
- ㉔ HDMI port
- ㉕ USB port
- ㉖ Connecting cable
- ㉗ Connection socket on endoscope
- ㉘ Connection socket on monitor

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Fragile

Keep dry

The device is subject to the requirements of the WEEE Directive 2012/19/EU.

Type BF device

Keep away from sunlight

Manufacturer

Humidity limitation

Atmospheric pressure limitation

Number of products in the product packaging

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Scope of protection against moisture:
Protection against continuous immersion

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1.2 Symbols for safety instructions

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

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2 Intended use

2.1 Intended use

The intubation videoscope with deflecting tip 10331 BX is used in combination with the C-MAC® monitor 8403 ZX and the C-MAC® PM 8403 XD for endotracheal intubation and inspection of the oropharynx.

Videoscopes are applied parts, which must only be used with KARL STORZ supply units which are intended for this purpose.

2.2 Contraindications

No contraindications relating directly to the product are currently known. The intubation videoscope must not be used if, in the opinion of the responsible physician, the patient's health is endangered through its use, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

2.3 Target groups

This instruction manual is intended for physicians, medical support staff and sterilization assistants.

2.4 User qualification

The intubation videoscope from KARL STORZ may only be used by persons with an appropriate medical qualification and who are acquainted with laryngoscopic techniques. The information given in these instructions only serves to instruct in the correct handling and preparation of the intubation videoscope. It is not intended as an introduction to the techniques of laryngoscopy.

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**2.5 Training in the operation and
function of the device**

Your local representative or responsible KARL STORZ member of staff is available to provide training and information on further training alternatives.

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**2.6 User profile of physician and
assistants**

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in one of the languages used on the device and in the instruction manual.
- Be thoroughly trained in the operation of the instruments.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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3 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

- Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

Combination of system components

- Combinations of medical devices are only assured to be technically safe if
 - they are listed as such in the respective instruction manuals or
 - the intended purpose and interface specifications of the devices used in combination permit this.The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.
- No high-frequency instruments may be used in combination with the intubation videoscope.
- Only KARL STORZ devices may be used. Disregarding this warning can endanger the patient and user.
- Only use accessories (tube holders, etc.) from KARL STORZ. Accessories from other manufacturers can damage the intubation videoscope.

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Risks of injury

- i** **NOTE:** All users must have been instructed in the operation and use of the C-MAC® VS, otherwise this endoscope may not be used.
- i** **NOTE:** The C-MAC® VS must be used in accordance with the airway algorithm in the respective country of application or the official specifications for airway management.
- Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.
- Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope. Always wear standardized protective equipment when using lasers.
- The high-intensity light emitted by the C-MAC® VS can result in increased temperatures at the light outlet. However, these always stay below the value which could cause burns.
- Risk of explosion if flammable anesthetics are used in the immediate proximity of the flexible intubation videoscope. The intubation videoscope must not be used in the presence of flammable or combustible media.
- If the flexible intubation videoscope malfunctions during use on the patient, discontinue use immediately. Place the distal tip of the flexible intubation videoscope in the neutral, non-deflected position and remove the intubation videoscope slowly and carefully from the patient.
- Danger of electric shock! Have servicing on these devices carried out by authorized personnel only. Any opening of the device by

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unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.

- Before performing any maintenance and cleaning work on the intubation videoscope, remove the connecting cable 8403 X and the tube holder.
- The intubation videoscope must not be used during the discharge of a cardiac defibrillator. Remove the intubation videoscope from the operation area before discharge takes place.
- Do not allow any live components or instruments to come into contact with the intubation videoscope.

Risks of infection

- The intubation videoscope is not sterile when delivered, and must therefore be properly cleaned and disinfected prior to initial use and each subsequent reuse.
- Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

Damage to the medical device

- Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.
- Intubation videoscopes must never be steam sterilized (autoclaved).
- The maximum permissible temperature for decontamination, sterilization, transport and storage is 65°C.

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- Also do not place the intubation videoscope on other devices that become very hot. Do not expose the device to direct sunlight or excessive heat, to avoid damage to the device.
- Intubation videoscopes must never be exposed to X-rays, to avoid damage.
- Do not hit the distal tip of the flexible intubation videoscope against hard objects, and never use a sharp object to remove dirt.
- Never lay heavy objects on the intubation videoscope.
- To avoid damaging the monitor 8403 ZX, do not remove the SD card or connecting cable whilst video recording.
- Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.
- The intubation videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for intubation videoscopes that are repaired by unauthorized personnel.

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4 Handling

4.1 Unpacking

WARNING: The intubation videoscope is not sterile when delivered, and must therefore be properly cleaned/disinfected prior to initial use and each subsequent reuse.

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

4.2 Basic equipment

The set 10331 BXK comprises:

- 1x C-MAC® VS 10331 BX
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025 E
- 1x Tube holder 10331 BAX

4.3 Possible combinations

The C-MAC® VS can be operated with the following monitors:

- C-MAC® motor 8403 ZX (software version: 704v300 or higher)
- C-MAC® PM 8403 XD
- Via the camera control unit C-HUB® II 20 2903 01 (software version: 721v105 or higher) with other monitors

The intubation videoscope must be used according to the medical rules and procedures recognized for laryngoscopy.

NOTE: All users must have been instructed in the operation and use of the C-MAC® VS, otherwise this endoscope may not be used.

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4.4 Preparation before use

NOTE: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding preliminary cleaning.

WARNING: Whenever the intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged intubation videoscopes or damaged accessories must not be used.

4.5 Inspecting the surfaces

WARNING: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

CAUTION: Do not squeeze the distal tip of the flexible intubation videoscope too tightly.

Every time the intubation videoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition. Intubation videoscopes with sharp edges must not be used. They may cause serious injury to the patient or user.

Check the connecting cable 8403X for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for loose connections.

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4.6 Holding the device correctly

CAUTION: The intubation videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the intubation videoscope.


Take the handle of the intubation videoscope in your right or left hand as shown in the image. With this kind of control grip, the intubation videoscope can be operated with one hand.

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4.7 Assembling the adaptor for tube fixation

If necessary, push the adaptor for tube fixation from distal onto the C-MAC® VS and fix in place by turning the knurled screw clockwise.

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
4.8 Initial operation

4.8.1 Connection to the C-MAC® monitor 8403 ZX

i **NOTE:** The instruction manual of the C-MAC® monitor must be observed! (Art. no. 96076008DER).

Plug cable into monitor

To switch on the device, insert the connecting cable of the C-MAC® VS into the corresponding connection socket 14 or 13 on the C-MAC® monitor 8403 ZX.

i **NOTE:** The monitor 8403 ZX makes it possible to connect two CMOS video units with Lemo connectors. All CMOS videoscopes from KARL STORZ can be connected to both connections. Video units from other manufacturers cannot be connected.

i **NOTE:** If two cameras are connected, the user must select the camera for which the images are to be displayed in the menu (refer to the instruction manual for the C-MAC® monitor 8403 ZX, Art. no. 96076008DER).

Plug cable into C-MAC® VS

Connect the plug of the connecting cable 8403 X to connection socket of the C-MAC® VS 10331 BX. Note the visible and tangible lug and push together both components until you feel resistance.

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Switch on monitor

Switch on the monitor by pressing the power switch ②. The green pilot lamp on the monitor lights up.

The C-MAC® VS is ready for immediate use, and the C-MAC® VS is immediately ready for use. Once the multifunction button on the C-MAC® VS lights up blue, videos and images can be recorded.

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the intubation videoscope out of the patient's body.

NOTE: The connecting cable can be removed during operation and connected to another videoscope.

NOTE: Make sure that the cable is fully inserted. There must not be a gap.

WARNING: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other units at any time.

WARNING: Optical radiation poses a risk of injury to eyes. Never look directly into the light output at the distal tip.

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4.8.2 Connection to the C-MAC® PM 8403 XD

NOTE: The C-MAC® VS 10331 BX can only be connected to the C-MAC® PM 8403 XD. The 8401 XD model is not compatible.

NOTE: The instruction manual for the C-MAC® PM must be observed (Art. no. 96076020DER).

Connect the C-MAC® PM to the socket on the C-MAC® VS. The orientation pin on the plug and socket facilitates connection. Note the visible and tangible lug and push together both components until you feel resistance.

! WARNING: Optical radiation poses a risk of injury to eyes. Never look directly into the light output at the distal tip.

NOTE: Make sure that the plug is fully inserted. There must not be a gap.

NOTE: The monitor can also be opened up when connecting both components.

The C-MAC® VS is ready for immediate use. As soon as the multifunction button on the C-MAC® VS lights up blue, a blue information bar will appear on the C-MAC® PM monitor for a few seconds. It is now possible to record images and videos.

NOTE: It is not necessary to perform a white balance.

NOTE: The image focus is adjusted automatically. Manual adjustment is not necessary.

NOTE: The C-MAC® PM can be removed during operation and connected to another laryngoscope.

NOTE: After 10 minutes, the display of the C-MAC® PM turns off automatically and must be turned back on again by closing it and reopening it or by disconnecting and reconnecting it (switching it off and on again).

NOTE: It may be more difficult to read the display in bright light, e.g., in direct sunlight.

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4.8.3 Connection to the C-HUB® II

NOTE: The C-MAC® VS can only be connected to the C-HUB® II 20 2903 01 (software version 721v105 or higher). Older models are not compatible.

NOTE: The instruction manual of the C-HUB® II must be observed (Art. no. 96206234DER)!

Connect the plug of the connecting cable 8403 X to connection socket 3 of the C-MAC® VS 10331 BX. Note the visible and tangible lug and push together both components until you feel resistance.

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4.8.4 Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

The C-MAC® VS is ready for immediate use. The multifunction button on the C-MAC® VS lights up blue at the same time.

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4.9 Assembling the endotracheal tube

Push the endotracheal tube from distal onto the intubation videoscope.

NOTE: When threading on the endotracheal tube, the distal tip of the intubation videoscope must be encased by the tube.

WARNING: Before the threaded on endotracheal tube is positioned in the patient via the vocal folds, the distal tip must be brought into the neutral position as otherwise the instrument could be damaged. To do this, release the deflection lever.

WARNING: Risk of injury to the patient: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

CAUTION: Damage to the instrument. When removing the endotracheal tube, the distal tip must be in the neutral position.

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4.10 Handling the deflection mechanism

- !** **NOTE:** The neutral position of the distal tip is straight, without any deflection whatsoever.
- !** **WARNING:** Risk of injury to users: Before each use, check the functionality of the deflection lever and the connected distal tip. Defective intubation videoscopes must not be used.
- !** **NOTE:** Test deflection without the endotracheal tube threaded on.
- !** **WARNING:** Risk of injury to the patient: When inserting the intubation videoscope into the patient, the distal tip must be straight and not deflected.
- !** **WARNING:** Risk of injury to the patient: Control the deflection up to the stop and avoid rapid or jerky movement so as not to injure the patient or damage the instrument. The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases.
- !** **WARNING:** Before the threaded on endotracheal tube is positioned in the patient via the vocal folds, the distal tip must be brought into the neutral position as otherwise the instrument could be damaged. To do this, release the deflection lever.

The distal tip can be deflected.

- $\geq 90^\circ$ without tube
- $\geq 60^\circ$ with endotracheal tube threaded on

- !** **CAUTION:** Damage to the instrument. When removing the endotracheal tube, the distal tip must be in the neutral position.

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4.11 Multifunction button 'BlueButton' on the C-MAC® VS

The intubation videoscope features a multifunction button ④.

Color coding meanings for the multifunction button 'BlueButton'

The multifunction button 'BlueButton' ④ on the C-MAC® VS can be used to activate single image capture and video recording. Documentation with the multifunction button on the C-MAC® VS can be performed on the C-MAC® PM 8403 XD and C-MAC® monitor 8403 ZX.

The colors of the multifunction button have the following meanings:

Blue: permanently illuminated. The recording function is enabled.

Green: temporarily illuminated when recording an image, or flashes during video recording.

i NOTE: The instruction manual for the C-MAC® PM must be observed (Art. no. 96076020DER).

i NOTE: The instruction manual of the C-MAC® monitor must be observed (Art. no. 96076008DER)!

Approx. 5 seconds after each time the C-MAC® PM is switched on, a blue bar is displayed at the top of the screen for a few seconds. The remaining memory capacity in minutes is shown on the left-hand side, while the number of the directory in which the saved data are stored is shown on the right and the current software version in the middle.

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
4.12 Single image capture

Single image capture is activated when the multifunction button ④ is pressed down once briefly. The current image on the screen is saved on the SD card in the connected monitor 8403 ZX / C-MAC® PM 8403 XD. During single image capture, the multifunction button ④ briefly lights up green and then switches back to blue.

4.13 Video recording

Video recording is started by pressing and holding the multifunction button ④ (for approx. 2 seconds) and is stopped by briefly pressing the multifunction button once more. The video recording is saved on the SD card in the connected C-MAC® monitor 8403 ZX or on the internal storage of the C-MAC® PM 8403.

During video recording, the multifunction button ④ on the C-MAC® VS flashes green, switching back to blue when pressed again briefly at the end of recording.

④ NOTE: The video stream can only be viewed on a PC if an MPEG 4-Codec is installed on the PC.

④ NOTE: If memory capacity of the C-MAC® PM 8403 XD is less than 10 minutes, the bar is displayed in orange-red.

4.14 Disassembly

Disconnect the connecting cable 8403 X from the C-MAC® monitor or C-HUB II.
Disconnect or separate the connecting cable or C-MAC® PM from the C-MAC® VS.
If necessary, release the adaptor for tube fixation from the C-MAC® VS by turning the knurled screw counter-clockwise. Only guide the tube holder cautiously over the deflection if the knurled screw is undone.

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- This reprocessing step is part of the validated process.
- Only for serial numbers from PR 1191: This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

NOTE: Empty fields are not part of the validated process.

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- !** **CAUTION:** Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.
- i** **NOTE:** National laws and regulations must be observed.
- i** **NOTE:** We recommend using a plastic basin for disinfection to avoid damaging the medical device. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.
- i** **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.
- i** **NOTE:** In principle, the C-MAC® VS and C-MAC® connecting cable 8403 X are resistant to manual and machine cleaning and disinfection. If the protective cap is not being used, then manual rinsing and drying of the electrical contacts may be necessary in order to ensure that functioning is not impaired. For this reason, the protective cap 8403 YZ is also offered. The user is responsible for deciding whether or not this cap needs to be used, depending on the individual process management (drying).
- i** **NOTE:** For information on reprocessing of the C-MAC® monitor 8403ZX, please refer to the corresponding instruction manual 96076008DER.
- i** **NOTE:** For information on reprocessing of the C-MAC® PM 8403 XD, please refer to the corresponding instruction manual 96076020DER.

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
Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Clean surface	27652
Ventilation connector	13242 XA
Contacts (optional)	8403 YZ

6 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

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- ⓘ **NOTE:** To avoid damaging the sheath, do not use pushing and pulling movements.
- ⓘ **NOTE:** Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

7 Leakage test

The medical device should be tested to ensure it is leakproof after use and before reprocessing.

7.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the medical device must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

- ⓘ **NOTE:** The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

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
Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the medical device is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

7.2 Performing the leakage test – in liquid

NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in the chapter 'Performing the leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid.

Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

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ATTENTION: In this case, fill out the accompanying repairs form and send the medical device, which has been cleaned by wiping to remove organic soiling/residue, for repair or, alternatively, contact KARL STORZ. The endoscope must be labeled as not disinfected. Line the storage/carrying case with the transport protection sheet 13990 SFN.

8 Disassembly

Remove the tube holder prior to cleaning and disinfection.

9 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility and necessary exposure time for the disinfectant must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

NOTE: Validation was performed with the Tristel Trio Wipes System from Tristel Solutions Limited.

NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

9.1 Manual reprocessing

9.1.1 Manual cleaning

The endoscope must be completely immersed in a cleaning solution. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush.

Finally, it must be rinsed with cold water to ensure neutralization.

ATTENTION: Never remove dirt with a needle or other sharp objects, since these can cause severe damage.

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ATTENTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

9.1.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

9.1.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

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9.2 Drying

Finally, the surfaces are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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10 Machine reprocessing

The following method for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

10.1 Connecting

In order to ensure effective machine cleaning and disinfection, the working channel of the endoscope must be connected up to the washer and disinfectant.

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10.2 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

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

11 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

i NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

i NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on the medical device, it must be dried again manually. The surfaces (e.g., connections) must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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- Afterwards, a functional check must be carried out.

i NOTE: Check for a clear view and steerability.

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13 Sterilization

WARNING: Before sterilization, all parts of the medical device must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

ATTENTION: The medical device must never be steam sterilized (autoclaved).

ATTENTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

-  **NOTE:** All deviations from the recommended sterilization parameters must be validated by the user.
-  **NOTE:** When using reprocessing equipment, please observe the equipment manufacturers' instructions.
-  **NOTE:** Sterilization is not possible on surfaces which have been lubricated and oiled.

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13.1 From SN PR 1191 onwards: Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

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The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device **from serial number PR 1191** onwards:

- STERRAD® 100S Long Cycle
- STERRAD® NX® Standard Cycle
- STERRAD® 100NX® Standard Cycle

13.2 From SN PR 1191 onwards: Hydrogen peroxide (H₂O₂) sterilization – STERIS® V-PRO™

NOTE: The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

The following STERIS V-PRO® sterilization methods have been validated and approved by KARL STORZ for this medical device **from serial number PR 1191** onwards:

- V-PRO® 1 Plus non-lumen Cycle
- V-PRO® 1 lumen Cycle
- V-PRO® 1 Plus lumen Cycle
- V-PRO® maX lumen Cycle
- V-PRO® 60 Flexible Cycle

13.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30-45 minutes.

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15 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes.

ATTENTION: Do not expose the medical device to direct sunlight or excessive heat.

ATTENTION: The distal tip of the medical device must be in the neutral position during storage.

ATTENTION: Do not store the medical device in the case.

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10331 BXH

15.1 Storage

- ! **CAUTION:** Do not expose intubation videoscopes to direct sunlight or excessive heat.
- ! **CAUTION:** The distal tip of the intubation videoscope must be in the neutral position during storage.
- ! **CAUTION:** The sheath of the intubation videoscope must never be kinked or coiled up too tightly. Intubation videoscopes should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage. It is most suitable to store the intubation videoscope suspended vertically with the distal end pointing down. The room where it is stored should be at a normal temperature and humidity.

The intubation endoscope can be safely attached to the mobile stand 8401 YA using the holder 10331 BXH (see Figure). It is secured directly to the cart's equipment rail.

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16 Reprocessing the tube holder

16.1 General warnings

WARNING: Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

WARNING: Risk of infection from non-sterile delivery. Reprocess the medical device before each use. Check the medical device for visible soiling before use.

WARNING: Before reusing the device, replace any defective parts or send the medical device for repair.

ATTENTION: Never use metal brushes to clean the medical device.

ATTENTION: Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

NOTE: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

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Accessories

Accessories required for carrying out reprocessing:

Accessories: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Clean surface	27652
Lumen	27650 A

17 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use.

Wipe-clean the exterior surfaces of the medical device with a damp cloth.

i NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

18 Disassembly

The tube holder must be removed from the endoscope prior to cleaning and disinfection.

19 Manual reprocessing

19.1 Manual cleaning

The medical device must be completely immersed in a cleaning solution. The lumina must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The lumina are brushed until the brush head is visibly clean. Never brush back and forth within the channel.

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19.1.1 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. The lumina must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

19.1.2 Manual High Level Disinfection (HLD)

The medical device must be completely immersed in a high level disinfectant solution. The lumina must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The medical device can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

19.2 Drying

Finally, the surfaces and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION: The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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20 Machine reprocessing

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

i NOTE: Both machine decontamination methods (chemothermal and thermal decontamination) are possible.

20.1 Connecting

In order to ensure effective machine cleaning and disinfection, the medical instrument must be connected up to the washer and disinfecter.

The selection of a suitable device holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off manually afterwards.

20.1.1 Machine cleaning / chemical-thermal disinfection

A maximum temperature of 65°C and a maximum pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

20.1.2 Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces and lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.
- Afterwards, a functional check must be carried out.

22 Packaging

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953).

23 Sterilization

! WARNING: Before sterilization, all parts of the medical device must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

! WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

i NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

i NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

i NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

23.1 Steam sterilization using the fractionated prevacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C – 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

or

- 134°C – 137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

23.2 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

ATTENTION: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

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The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle
- STERRAD® NX® Standard Cycle
- STERRAD® 100NX® Standard Cycle

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The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® 1 Plus non-lumen Cycle
- V-PRO® 1 lumen Cycle
- V-PRO® 1 Plus lumen Cycle
- V-PRO® maX lumen Cycle
- V-PRO® 60 Flexible Cycle

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23.5 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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24 Technical description

24.1 System description

The intubation videoscope from KARL STORZ has a working length of 410 mm and a diameter of 5.5 mm, and is made of stainless steel and a flexible Pellethane sheath. The PEEK (polyetheretherketone) distal tip is 5.5 mm in diameter and houses the CMOS video module. The distal tip of the intubation videoscope can be deflected through $\geq 90^\circ$ without the tube. The handle houses the deflection lever and the function button.

The intubation videoscope from KARL STORZ is connected to the recommended C-MAC® monitor 8403 ZX. If an SD memory card is inserted in the C-MAC® monitor 8403 ZX, images and videos can be saved using the function button on the videoscope and the monitor which are assigned the functions of video recording and image capture.

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24.2 Technical data

24.2.1 C-MAC VS 10331 BX

Direction of view:	0°
Angle of view:	100°
Deflection:	≥ 90° without tube ≥ 60° with tube

Working length:	41 cm
Outer diameter:	5.5 mm
Illumination:	Integrated LED light source

Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight	445 g
Drip water protection:	IPX8
Type of protection against electric shocks:	Protection Class II

Provides Type BF protection when used in conjunction with a KARL STORZ camera control unit with a BF symbol on the connection port.

Operating/storage conditions

Operating conditions:	0°C to +40°C 30% to 70% rel. humidity
Storage conditions:	-20°C to +50°C 5% to 95% rel. humidity
Atmospheric pressure:	700 to 1080 hPa

CAUTION: The specifications in relation to atmospheric pressure can only be guaranteed if the pressure equalization cap is fitted to the vent port during transportation.

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25 Service, repair and disposal

25.1 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

25.2 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

25.3 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

25.4 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see Chapter 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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25.5 Electromagnetic Compatibility

The applied part described in this instruction manual has been tested as a system with the following basic devices. For relevant electromagnetic compatibility (EMC) information, refer to 'Electromagnetic Compatibility (EMC) information' in the instruction manuals for these basic devices

Basic device	Instruction manual
C-MAC® monitor 8403 ZX	96076008 DER
C-MAC® PM 8403 XD	96076020 DER
C-HUB® II 20 2903 20	96206529 DER

The EMC warning statements, precautions, notes and Emission/Immunity limits specified in the instruction manuals for the basic devices also apply to the applied part described in this instruction manual.

25.6 Directive compliance

In accordance with the Medical Device Directive (MDD):
C-MAC® VS 10331 BX:

This medical device belongs to Class I

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

25.7 Important information on dispatch

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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25.8 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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26 Appendix

26.1 Accessories

8401 YA: Stand,
height 120 cm, rollable, with five legs and anti-static castors, crossbar 25 cm x diameter 25 mm, for positioning the monitor, with tray for laryngoscopes, dimensions (w x d x h): 30 x 20 x 10 cm

8401 YAA: Crossbar,
for stands 8401 YA, 50 cm x diameter 25 mm, for positioning the C-MAC® monitor 8403 ZX

8401 YAB: Crossbar,
for stands 8401 YA, 70 cm x diameter 25 mm, for positioning the C-MAC® monitor 8401 ZX

8401 YB: Multifunctional clamp,
with VESA 75 standard mount. Suitable for mounting on tubes up to 24 mm in diameter and rails up to 25 mm in height

10331 BXH: Holder for C-MAC® VS,
endoscope holder for C-MAC® rigid intubation videoscope 10331 BX, made of Plexiglas, distal open end, for securing to a standard rail

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
200200 85: Airway TROLL-E mobile stand,
on 4 antistatic double castors, 2 of which
are lockable, with VESA 75/100 holder for
mounting a monitor, integrated cable channel
in the central upright, max. load capacity of
monitor module 15 kg, dimensions of monitor
cart: 670 x 1660 x 670 mm (W x H x D), castor
diameter: 100 mm.
Mobile stand is not assembled on delivery.

200200 86: C-MAC® TROLL-E,
on 4 antistatic double castors, 2 of which are
lockable, with stainless steel tube,
dimensions of mobile stand: 670 x 1500 x
670 mm (w x h x d),
castor diameter: 100 mm
Mobile stand is not assembled on delivery.

29005AK: Storage basket,
for accessories, for mounting to equipment rail
29003 GS,
dimensions (w x h x d): 300 x 100 x 200 mm,
for use with equipment cart 29005

20020047: Console,
incl. installation accessories, max. capacity
12 kg, dimensions: 455 x 350 mm (w x d), for
use with subracks for mobile Stands 20020060
and 20020061

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Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +84 28 3823 8000
Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Janguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 56381888, Fax: +86 10 56381999
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 1595, Australia
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)

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