TPO

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USER'S MANUAL

SelfyCheck Prim Test Strips

produced by Taidoc Technology Corporation

Before use

Dear User of SelfyCheck Prim Test Strips!

Read this User's Manual carefully before using the product.

Pay special attention to all warnings and precautions that are marked with  sign.

All elements of the System are designed, tested and checked for compatibility in such a way as to ensure accurate results of measurement of blood glucose level.

 *Use SelfyCheck Prim Test Strips only as a part of the system with SelfyCheck Prim Glucose Meter and SelfyCheck Prim Control Solutions.*

Regular glycemic control allows you and your doctor to monitor the diabetes in a better way. Due to the small size and simplicity of operation the System provides a possibility of self-monitoring of the blood glucose level anywhere and at any time.

Many people, performing the blood analysis themselves, find it useful to start with practising on the control solution (see Section "CHECK OF SYSTEM OPERATION USING THE CONTROL SOLUTION").

Please contact the support service or place of purchase regarding the issues related to the use of this device.

1. MEDICAL DEVICE NAME

SelfyCheck Prim Test Strips (the "Device", "Medical Device", "Test Strips").

2. MEDICAL DEVICE COMPOSITION

Medical Device **SelfyCheck Prim Test Strips** includes the following versions:

Version 1:

1. SelfyCheck Prim Test Strips (10 pcs. in a pack) – 1 pce.
2. User's Manual – 1 pce.

Version 2:

1. SelfyCheck Prim Test Strips (25 pcs. in a pack) – 1 pce.
2. User's Manual – 1 pce.

Version 3:

1. SelfyCheck Prim Test Strips (50 pcs. in a pack) – 1 pce.
2. User's Manual – 1 pce.

Version 4:

1. SelfyCheck Prim Test Strips (100 pcs. in a pack) – 1 pce.
2. User's Manual – 1 pce.

3. INFORMATION ON THE MANUFACTURER AND AUTHORIZED REPRESENTATIVE OF MEDICAL DEVICE

Manufacturer: Taidoc Technology Corporation

Manufacturer's address: B1-7F, No.127, Wugong 2nd Rd., Wugu Dist., New Taipei City 24888, Taiwan, R.O.C. Tel. +886-2-6625-8188, E-mail: sales@taidoc.com.

Manufacturing site: 1F & 2F & 3F, No. 12, Lane 5, Sec. 2, Nanshan Rd., Lujhu Township, 33852 Taoyuan County, Taiwan.

Authorized Representative (complaint management organization):

NEXT Bio LLC

Address of the authorized representative: 5 Akademika Pavlova Street, lit. B, room 46-H, 197022 Saint Petersburg, Russia

Tel.: +7 (812) 703-79-75 (multiline), fax: +7 (812) 703-79-76

Hot Line: 8-800-333-43-76 (free throughout Russia)

4. PURPOSE

The test strip is intended to quantify the level of glucose for use as a support for monitoring glucose levels in fresh capillary whole blood (taken from the finger, palm, forearm, upper arm) as a part of the SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter) for self-testing by people with diabetes at home and at stay in clinical settings.

5. TESTING PRINCIPLE

The device is intended for individual use.

SelfyCheck Prim Test Strips are used as a part of SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter). For check of system operation it is necessary to use SelfyCheck Prim Control Solutions with low, normal and high concentration from the starting set included in device components.

User's Manual is supplied together with the device to ensure correct use of the device, self-study and learning of skill of self-monitoring of blood glucose level.

Testing Principle

Determination of glucose in the blood sample is based on the measurement of electric current induced by the reaction of glucose and enzyme on SelfyCheck Prim Test Strip.

After inserting the test strip into the blood glucose meter test port the System is brought to the blood sample with the free end of the test strip. The blood is automatically drawn into the absorption channel where the reaction takes place. The transient electric current is generated during reaction. The system reads the current strength value obtained during the reaction between the enzyme and glucose in the blood sample in the reaction cell of the test strip and calculates the concentration of blood glucose level using the pre-set equation of the software built into the blood glucose meter. The current strength value depends on the blood glucose concentration, batch of test strips, operating voltage and temperature. The result is shown on display.

Sample of fresh capillary whole blood for measurement of blood glucose level can be obtained not only from the fingertip, but also from other body parts.

6. INDICATIONS

Main indications for usage of the medical device – diabetes mellitus type 1 and 2.

7. CONTRAINDICATIONS

There are no absolute contraindications. There are use limitations.

8. POSSIBLE SIDE EFFECTS

Not identified.

9. SCOPE

Medical Device “SelfyCheck Prim Test Strips” are used in the endocrinology.

10. POTENTIAL USERS

The device can be used for self-monitoring by people at home as well as by health care professionals for control of effectiveness of measures for self-monitoring of diabetes in patients of medical and preventive treatment facilities.

The device can be used by lay users with various physical and mental abilities not having special medical, scientific, and technical knowledge with regard to blood glucose level monitoring, having read the User’s Manual in advance.

11. DEVICE USE LIMITATIONS

- ⚠ The Test Strips are designed for external use only (for in vitro diagnostics).
- ⚠ Test strips as a part of the System are not recommended to use in dry environment, especially with synthetic materials (synthetic clothes, carpets, etc.). It can lead to the damage of the System by static discharge and, therefore, to the erroneous results.
- ⚠ It is prohibited to use the test strips as a part of the System in close proximity to the sources of strong electromagnetic radiation.
- ⚠ If the System is used at home, the measurement is performed only on the samples of fresh capillary whole blood; in the medical institutions the healthcare professionals can use the samples of fresh capillary and venous whole blood.
- ⚠ The device is not intended for use by persons with special needs on their own (children under 18, disabled persons, the elderly). Their usage of the device requires assistance and help.
- ⚠ The device shall not be applied to patients with severe infections, malignant neoplasms, massive oedema.
- ⚠ The device is not intended for diagnosis of diabetes mellitus.
- ⚠ The device is not designed for measurement of blood glucose level in the new-born infants.
- ⚠ Possible factors of impact on the results of analysis:
 - ⚠ The checked blood was squeezed out of finger by strong press and therefore contains the interstitial fluid.
 - ⚠ Shelf life of test strips has expired (shelf life is 21 months from the date of manufacture or 6 months from the date of container opening).
 - ⚠ Test strips were stored at improper temperature (proper storage temperature: from +2°C to +30°C, not more than 85% of relative humidity).
 - ⚠ Insufficient moisture protection of test strips.
 - ⚠ Sanitizer, used in the preliminary disinfection of the point of blood sampling, not completely vaporized from the skin.
 - ⚠ The measurement was performed with unwashed hands.
 - ⚠ After washing the hands were not dry for further measurement of blood sugar level.
 - ⚠ The excessively high blood pressure can lead to the distortion and erroneously low measurement results.

- ⚠ Very high and very low hematocrit (lower than 35% or higher than 60%) can lead to the incorrect results. Hematocrit value shall be within the range from 35% to 60%.
- ⚠ The results can be inaccurate in case of the patients with very low arterial blood pressure or those who are in shock.
- ⚠ Endogenous substances, naturally present in the human body (bilirubin, uric acid, glutathione, triglycerides, other medicamental chemical compounds), as well as medicinal products do not have a significant impact on the results of blood glucose level measurement using a blood glucose meter if normal concentrations are not exceeded or in case of standard therapeutic doses. High concentration in blood can affect the analysis results. It is recommended to pay special attention to the interpretation of these results.
- ⚠ In case of patients with impaired peripheral blood circulation, it is not recommended to collect the blood from alternative body parts, because the obtained readings can differ from actual blood glucose level. Such non-conformities can be observed in the following cases: In case of severe dehydration as a result of diabetic ketoacidosis, hyperglycaemia due to the stress, hyperosmolar nonketotic com, shock, decompensated heart failure of class IV according to NYHA scale or peripheral arterial occlusive disease. Measurement errors can lead to the incorrect interpretations and thus to the serious health consequences.

12. APPLICATION CONDITIONS

The device is intended for application in hospitals, outpatient departments, medical institutions, ambulance cars as well as at home and in the field.

13. APPLICATION PRECAUTIONS

- ⚠ The test strips shall be supplied and stored in the plastic container with tightly closed cover. The container cover shall be tightly closed immediately after taking out one test strip.
- ⚠ The test strips are for single-use, their re-use is forbidden and may cause cross infection.
- ⚠ Antiseptic rules shall be followed when using the device. It is necessary to wash hands with soap in warm water or use wet disinfecting wipe and let the skin dry before use.
- ⚠ It is necessary to follow standard precautions when handling blood. All samples and materials that came in contact with the patient are considered to be biohazardous, they shall be handled as potential infection carriers.
- ⚠ The device contains small details that may cause asphyxia and injuries in young children in case of inhalation and swallowing. The device shall be kept out of reach of children!
- ⚠ Do not remove the test strip from the test slot of the blood glucose meter during the analysis.
- ⚠ Self-measurement of blood glucose level does not replace the diagnosis by the doctor or regular appointments. Self-testing of blood glucose level shall be used only after receiving the detailed instructions for monitoring and treatment of diabetes from the qualified medical specialist (endocrinologist).

14. DESCRIPTION AND SPECIFICATIONS

14.1 SelfyCheck Prim Test Strips

⚠ The test strip consists of the following parts:

1. Absorbent Hole

Area of contact of test strip and blood sample. In case of contact the blood is absorbed automatically.

2. Confirmation Window

The filling of the window confirms that there is enough blood for analysis.

3. Test Strip Handle

Area for holding the test strip when inserting into the test slot.

4. Contact Bars

Area of contact of the test slot and test strip. When inserting the test strip into the test slot, the blood glucose meter will switch on automatically.

⚠ The test strip shall be inserted firmly until it will go no further with front side facing up.

⚠ If contact bars are not fully inserted into the blood glucose meter, the test result can be incorrect.

⚠ The test strips shall be used within 6 months from the date of first opening of the vial (but not later than 21 months from the date of manufacture). It is necessary to record the date of first opening of the container on its label.

⚠ The container shall be tightly closed immediately after taking out a test strip.

⚠ Touch the test strip only with clean and dry hands. Use each test strip immediately after taking it out of the container.

⚠ Do not use the test strip after expiration of the shelf life. It can lead to the incorrect results.

⚠ Do not bend, cut or alter the test strip.

Specifications


Model	TD-4311
Chemical components	Enzyme - GDH-FAD 8% Potassium phosphate 8% Ferricyanide 55% Carboxymethyl cellulose 29%
Dimensions of test strip:	7.5 x 36.06 x 0.5 (±0.2) mm
Dimensions of container:	53.0(±0.5) mm, Ø 30.46(±0.25) mm.

Weight of one test strip:	0.17 (± 0.003)g.	
Gross weight of container:	- 10 pcs. in container:	17.8 g $\pm 10\%$
	- 25 pcs. in container:	20.4 g $\pm 10\%$
	- 50 pcs. in container:	24.8 g $\pm 10\%$
Sample size:	0.5 μ l	
Measurement units:	mg/dl or mmol/l	
Reaction time:	5 s	
Measurement range:	from 1.1 to 36.1 mmol/l (from 20 to 650 mg/dl)	
Analyte:	Dextroglucose	
Enzyme:	GDH-FAD enzyme (glucose dehydrogenase – flavin-adenine dinucleotide)	
Hematocrit range:	35% – 60%	
Shelf life:	21 months from the date of manufacture or 6 months from the date of container opening.	
Sample type:	fresh capillary whole blood; whole dark blood	

14.2 Brief Description of Devices Intended for Use in Combination with the Stated Medical Devices

14.2.1. SelfyCheck Prim Blood Glucose Meter

Figure

Model:	TD-4183
Power source:	1.5V AAA alkaline battery
Battery life:	at least 500 measurement cycles
Display:	LCD, white illumination
Memory:	1,000 measurement results with date and time
External output:	Cable (strip port)
Enzyme type	GDH-FAD enzyme (glucose dehydrogenase – flavin-adenine dinucleotide)
Analysis time:	5 s
Measurement units:	mg/dl or mmol/l
Measurement range:	20–650 mg/dl (1.1–36.1 mmol/l)
Resolution (measurement interval):	0.1 mmol/l (1 mg/dl)

Sample size:	0.5 µl
Hematocrit level:	35–60 %
Consumed power:	157.08 mW
Dimensions:	52.3 x 90.3 x 18.0 (±0.3) mm
Weight (with battery):	57.4 g ± 15%
Operating mode setting time:	90 sec
Software Version	2011_0309 PM 02:08
Additional functions:	Automatic detection of sample Automatic countdown start Automatic shutoff after 3 minutes with no activity Warning of violation of the temperature conditions

14.2.2 SelfyCheck Prim Control Solutions

SelfyCheck Prim Control Solution is intended for usage in combination with SelfyCheck Prim Test Strips. SelfyCheck Prim Control Solution is used for calibration of the SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter). The device is intended for in vitro diagnostics only.

Model	SelfyCheck Prim (Low) Control Solution			TD-4907
	SelfyCheck Prim (Norm) Control Solution			TD-4908
	SelfyCheck Prim (High) Control Solution			TD-4909
Nominal vial volume	4 ml			
Vial dimensions	58.2(±0.6) mm Ø 17.5 (±0.35) mm			
Weight	6.6 g ± 10%			
Composition	<div>Level</div> <div>Composition</div>	Low	Normal	High
	Dextroglucose	0.05% (2.9 ±0.5 mmol/l)	0.12% (6.7 ±1.0 mmol/l)	0.25% (14.7 ±2.2 mmol/l)
	Normal saline	99.9%	99.83%	99.7%
	Amaranth	0.05%	0.05%	0.05%

Do not use control solution after the expiration date or 3 months after the initial opening of the vial. Note the date of opening on the vial with the control solution and dispose of the remaining solution after 3 months.

It is recommended to perform the control solution test at the temperature from 20°C to 25 C and relative humidity from 30% to 85%. Make sure that the control solution, blood glucose meter and test strips are within the specified temperature range before testing.

Shake the vial before use. To avoid bubbles do not shake the vials too vigorously.

15. OPERATION

15.1 Check of the Proper Combined Operation of the Blood Glucose Meter and Test Strips Using SelfyCheck Prim Control Solutions

It is recommended to perform the independent check of the joint operation of the blood glucose meter and test strip using SelfyCheck Prim Control Solution (supplied in the stating set as a part of the System or purchased separately). The control solution is put into the vial with the seal cap. The vial with "Low" marking is required for check and display of the low glucose level; the vial with "Norm" marking is required for check and display of normal glucose level; the vial with "High" marking is required for check and display of high glucose level.

The allowable range of results in the control solution analysis is indicated on the container with test strips. To obtain accurate results, it is necessary to conduct such diagnostics regularly.

The control solution test shall be performed:

- if you suspect that the blood glucose meter or test strips work improperly;
- if your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate;
- for practising the testing process;
- if you have dropped or think you may have damaged the device.

Performing a Control Solution Test

For performance of a control solution test you will need the blood glucose meter, test strips and control solution.

1. Insert a test strip

Insert the test strip into the slot with the contact strips forwards and upwards. The Blood Glucose Meter will turn on automatically with "☰" and "🔋" symbols appearing on the display.

2. Mark the test with the control solution

Press the "M" button to mark the test with the control solution. When the "QC" symbol is displayed on the device, the test results are stored in memory. If you decide not to carry out the test with the control solution, press the "M" button again and the "QC" symbol will disappear.

 *When conducting a test with a control solution, it is necessary to switch to the "QC" mode so that the result is NOT stored in the device memory along with the TEST RESULTS with blood. Otherwise, the results of the blood glucose test will be mixed in the device memory with the results of the control solution test.*

3. Bring the test strip to the drop of control solution

Turn the closed vial over several times (do not shake) before opening. The drop of the control solution shall be applied to the clean vial cap; then the test strip, inserted into the blood glucose meter, shall be brought to the drop. Once the absorbent hole of the test strip touches the drop of the control solution, the drop will be drawn in automatically.

! To avoid contamination of the control solution, do not apply the solution directly to the test strip.

4. Read and verify the result

The confirmation window of the test strip will be filled in case of sufficient sample size; after that the blood glucose meter will produce a sound signal and start the test. After the end of the countdown from 5 to 0, the test result with the control solution will appear on the display. Compare the result with the range indicated on the test strip container. The result must be within this range.

If the result is outside the specified range

If the repeated test results continue to fall outside the range specified on the container, the device can work or be used improperly, or the blood glucose meter and test strips can be incompatible. **DO NOT USE** this kit for testing of a blood sample. Contact the support service or place of purchase for help.

! The range of values for the test with the control solution, indicated on the container of test strips, is intended only for evaluating the results of the test with the control solution. This range is not a recommendation for blood glucose.

15.2 Performance of Blood Sample Test

For performance of the blood sample test you will need:

1. SelfyCheck Prim Blood Glucose Meter
2. SelfyCheck Prim Test Strips
3. SelfyCheck Prim Control Solutions (Low, Norm, High)
4. Lancing device (with AST tip).
5. SelfyCheck Lancets.

! Complete set of the System is specified on the package.

Components that can be purchased additionally:

1. SelfyCheck Prim Test Strips


To reduce the risk of infection:

! never share a lancet or the lancing device with other people;

! always use new sterile lancet. Lancets are for single use only;

! avoid getting lotions, oils, dirt and debris inside or on the lancets and the lancing device.

1. Setting up the Lancing Device for Blood Testing


1) Turn the lancing device cap and remove it by pulling it up.

2) Insert the lancet into the lancet holder and press it in order to fix it firmly.


3) Remove the protective cap from the lancet.


4) Put the device cap back and turn it until clicking sound to fix it (see figure).

The cap is fixed if the fixing indicator is in the extreme right position.

5) Hold the cap and choose the puncture depth turning the adjustable tip in any direction so that the arrow on the cap points at the required puncture depth.


Hold the device by the housing and not by the cap.

6) Pull by the slope adjustment area until click. When the device is ready for use the release button colour will be changed.

 If you did not hear the click possibly the device jammed before the lancet fixation.

2. Insert the Test Strip into the Blood Glucose Meter to Switch It On

Wait for “” and “” symbols to appear on the display.

3. Select the Suitable Measurement Mode

The device has four measurement modes: General, AC, PC and QC.


MODE	CONDITIONS FOR USE
General	any time of day, without regard to timing of food ingestion
AC	8 hours before meals minimum
PC	2 hours after a meal
QC	diagnostics using the control solution

The switching between the modes is performed in the following way:


1. “▲” symbol or “Gen” will appear on the display after switching on.

2. Press “▲” button to switch the modes General, AC and PC.

3. Press “M” button to switch to QC mode.

4. Obtaining the Blood Sample

Preparation of the Puncture Site

Stimulating blood circulation has a significant impact on the quality of the obtained result of blood glucose level measurement; to do so rub the puncture site before blood draw. The value of the glucose level in the blood collected from the site that has not been rubbed preliminarily differs essentially from the value of glucose level in the blood collected from the finger. Rubbing the puncture site before the blood draw allows for significant reduction of this difference.

Please follow the instructions below before collecting a drop of blood:

- ⚠ wash and dry your hands before starting the procedure;
- ⚠ select the puncture site at fingertip or on other part of the body;
- ⚠ rub the puncture site for about 20 seconds before penetration;
- ⚠ clean the puncture site using disinfecting wipe and let the skin dry;

5. Blood Drawing from the Finger

- 1) Puncture the finger

Press the lancing device tip to the bottom part of the finger pad tightly.

In order to make puncture press the release button. After the puncture completion the click will sound.

- 2) Remove the device from the puncture site. Use the lancing device to prick the skin. After the puncture completion remove the first blood drop using clean cotton. It is necessary to obtain not less than 0.5 µl of blood sample. Press to the puncture site

carefully to collect another blood drop. Be careful in order NOT to smear the blood sample.

⚠ Choose a different spot each time you perform a test. Repeated punctures at the same spot may cause pain feelings and calluses.

⚠ Consult the doctor before collecting the blood from the alternative sites.

⚠ It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

5. Apply the Blood Sample

It is necessary to bring the test strip inserted in the glucose meter to the drop of blood at an angle. Once the absorbent hole of the test strip touches the drop of blood, the blood will be drawn in automatically. The confirmation window of the test strip will be filled in case of sufficient blood sample size; after that the blood glucose meter will produce a sound signal and start the test.

⚠ DO NOT REMOVE the finger until you hear the sound signal.

⚠ Do not press the puncture site to the test strip and do not smear the blood.

⚠ If you do not apply the blood sample to the test strip within 3 minutes, the device will switch off automatically. To start a new measurement remove and replace the test strip.

⚠ The blood shall completely fill the confirmation window to start the count down. DO NOT try to add more blood to the test strip after the first drop is absorbed. Remove the used test strip and repeat the measurement using new test strip.

⚠ If you have trouble filling the confirmation window, consult your doctor or support service.

6. Reading of Result

The result of blood glucose test will appear after the countdown from 5 to 0 and will be saved to the device memory automatically (100 mg/dl = 5.5 mmol/l).

⚠ *Severe dehydration and excessive water loss can lead to the readings lower than the actual values. Immediately contact your doctor in case of suspicion of severe dehydration.*

⚠ *If the blood glucose level readings are lower or higher than the normal value, and the corresponding symptoms of the disease are absent, repeat the test first. Follow the recommendations of your doctor in case of symptoms or if the repeated tests continue to show the results lower or higher than the normal value.*

⚠ *To measure the blood glucose level use only fresh whole blood samples. In case of use of other samples the results will be incorrect.*

⚠ *If you have the symptoms that do not correspond to the test results and subject to adherence to the instructions provided in this Manual.*

⚠ *Check, if the result of analysis is within the acceptable range of results. Refer to your doctor, if the value exceeds 14 mmol/l.*

The System allows obtaining a result equivalent to the value for blood plasma.

Time of measurement	Range of normal values of blood plasma glucose level for people without diabetes (mg/dl)
Fasting and before meal	< 100 mg/dl (5.6 mmol/l)
2 hours after a meal	< 140 mg/dl (7.8 mmol/l)

⚠ *Consult your doctor to define the target range of values appropriate for you.*

Blood or plasma glucose concentration level can be expressed in units of weight (mg/dl) or molar concentration (mmol/l). Principle of approximate conversion of mg/dl value to mmol/l value is as follows:

mg/dl	divide by 18	= mmol/l
mmol/l	multiply by 18	= mg/dl

For example,

1. $120 \text{ mg/dl} \div 18 = 6.6 \text{ mmol/l}$;


2. $7.2 \text{ mmol/l} \times 18 = 129 \text{ mg/dl}$ (conversion results are approximate).

7. Disposal of Used Materials

Removal of the used test strip

Remove the test strip by moving the glucose meter ejector pin slider to the end position; the glucose meter will switch off automatically.

Removal of the used lancet

1) Twist off the device cap and remove the lancet. Dispose the used lancet following precautions: put the lancet protective cap on a hard surface and insert the used needle in it.

2. In order to remove the used lancet move the ejector forward. Dispose the used lancet according to the regulations of local legislation.

⚠ *Used lancets and test strips are considered to be biohazard waste. Dispose the used lancet according to the regulations of local legislation.*


⚠ *Use transparent tip for blood draw from the alternative sites (not from the finger pad).*

15.4 Blood Draw from the Alternative Site

It is possible to carry out analysis of a blood sample drawn from the alternative sites.

⚠ *There are limitations in terms of possibility of blood draw from the alternative sites. Please consult your physician before using this method.*

Replace the cap of the delivery device with a transparent cap for blood draw from the alternative sites. Pull the cocking control back until it engages. When making injections to the forearm, upper arm or hand, avoid sites with well defines veins, to avoid excessive bleeding.

When should the method of blood draw from the alternative sites be used?

Dietary habits, treatment, physical activity level or stress and diseases may influence blood glucose level. Capillary blood in finger pads reflects such changes faster than blood in other parts of the body. Therefore, **only the blood samples from finger** should be used when measuring blood glucose level during or immediately after a meal, exercise or any other actions or events.

It is strongly recommended to perform blood draw from the alternative sites **ONLY** under the following conditions:


- Before a meal or on an empty stomach (not earlier than 2 hours after the last meal).
- Not earlier than 2 hours after the insulin administration.
- Not earlier than 2 hours after physical exercise.

DO NOT USE blood draw from the alternative sites in case:

- You assume low level of blood glucose.
- You are not sure whether you have hypoglycaemia.
- You are performing test for hyperglycaemia.
- The result of blood glucose level measurement from the alternative sites does not correspond to your state of health.
- The results of your ordinary checks are subject to fluctuations.

Alternative sites for the blood draw:

1. finger,
2. palm
3. forearm
4. upper arm .

16. DEVICE CONTACT WITH HUMAN BODY

The test strips contact with damaged skin for a short period of time.

When producing the device the components of animal origin are not used. The product does not contain substances that can be considered as drugs or pharmaceuticals. Materials of biological origin intended for recycling and re-use are not used in the device.

17. LABELING

Label on the container body with test strips contains the following information: name, model, trademark, name and address of the manufacturer, name and address of the authorised representative, hot line, serial number, place for indication of the date of opening, number of test strips in the pack, "Warning! Refer to the instructions for use" inscription; "Best before" inscription (in YYYY-MM format); "In vitro diagnostics medical device" inscription; "Temperature range" inscription; "Humidity range" inscription; "Re-use is forbidden" inscription; range for qualitative interpretation of measured values in mmol/l and mg/dl units (low, normal, high levels) for control solution (in tabular form); "Use only with SelfyCheck Prim Blood Meter and SelfyCheck Prim Control Solutions" inscription; "Use within 6 months after the initial opening" inscription; "Do not freeze" inscription; internal label code; date and number of the Marketing Authorisation, opening control (strip, glued to container for strips cover)-depending on the release form.

Test strips are produced in two forms of release:

- 1) SelfyCheck Prim Test Strips - with opening control
- 2) SelfyCheck Prim Test Strips - without opening control

Labels of the consumer packaging of the test strips contain the following information: name, model, trademark, name and address of the manufacturer, name and address of the authorised representative, hot line, serial number, place for indication of the date of opening, number of test strips in the pack, "Warning! Refer to the instructions for use" inscription; "Best before" inscription (in YYYY-MM format); "In vitro diagnostics medical device" inscription; "Protect from direct sunlight" inscription; "Protect from moisture" inscription; "Temperature range" inscription; "Re-use is forbidden" inscription; period of measurement to obtain the result (5 sec); "Use only with SelfyCheck Prim Blood Meter" inscription; "Use within 6 months after the initial opening" inscription; internal label code; hot line; date and number of the Marketing Authorisation.

Symbols used in the product labelling

Labeling symbol	Interpretation
	Manufacturer
	Device batch number according to the manufacturer's numbering system
 YYYY-MM	Use before the specified date
	Re-use is forbidden
	Storage temperature range
	Protect from moisture

	Refer to the User's Manual
	Warning! Refer to the User's Manual
IVD	In vitro diagnostics medical device

18. REQUIREMENTS FOR ENVIRONMENT PROTECTION

The medical device is safe and does not affect the environment. Potential biohazard: doctors or people who use this product for analysis for several patients should be aware that any items that come in contact with blood, even after treatment, can be sources of vector-borne viral infections. The product does not generate its own energy and does not absorb the user's energy. The product is not toxic at room temperature, does not release toxic substances into the environment, does not have any harmful effect on the human body through direct contact. The product does not contain liquid and gaseous components in its composition.

19. SAFETY INFORMATION

The manufacturer does not guarantee the correct operation of SelfyCheck Prim Test Strips when used with any glucose meter other than SelfyCheck Prim Blood Glucose Meter. Usage of other glucose meters together with SelfyCheck Prim Test Strips may lead to erroneous results. The product cannot be used until it has been tested with the SelfyCheck Prim Control Solution and the results obtained with the test are within the required range. To avoid cross-contamination, it is not recommended to use the same lancing device by several people. It is prohibited to use the lancet that has been used by another person. To avoid infection, a new sterile lancet is to be used every time. Any items that come in contact with blood, even after treatment, can be sources of vector-borne viral infections.

20. QUALITY STANDARDS

The device corresponds to the international standards: IEC 62366:2007, BS EN 62366:2008, ISO 14971, IEC 62366-1:2015, EN ISO 15197:2015, ISO 23640: 2011, EN ISO 17511, EN 13532:2002, EN ISO 15223-1:2014, EN ISO 18113-4:15, EN ISO 18113-5, EN ISO 13485, CE MDD (93/42/EEC), IEC 61326-1.

21. STERILITY

The medical device is supplied non-sterile.

22. METHODS AND MEANS OF DISINFECTION AND PRE-STERILISATION CLEANING

This device is intended for single use, cleaning for re-use is not acceptable.

23. STORAGE AND TRANSPORTATION CONDITIONS, SHELF LIFE

23.1 Transportation

Transportation can be carried out using all types of transport in covered vehicles or general-purpose containers protected from sunlight in accordance with cargo transportation rules applicable to them. Transportation conditions: within the temperature range of 2°C to 30°C and relative humidity up to 85% in the original container.

23.2 Shelf Life

Shelf life of the test strips in a closed container before the first opening constitutes 21 months from the date of manufacture. Shelf life of the test strips after opening of the container constitutes 6 months, but not more than 21 months from the moment of manufacture. Do not use the test strip after expiration of the shelf life.

23.3 Storage

Store the device at 2°C – 30°C and relative humidity up to 85%.

- Store the device in original container only, do not transfer to another package.
- In a dry cool place.
- Avoid direct sunlight and heat.
- Do not freeze.
- Keep the container with test strips away from children.

23.4 Use

Device is to be used within the temperature range from 10°C to 40°C and relative humidity not more than 85 %.

- The container cover shall be tightly closed immediately after taking out of one test strip.
- Touch the test strip only with clean and dry hands. Use each test strip immediately after taking it out of the container.
- Never bend, cut or damage the test strip.
- Do not use the test strip after expiration of the shelf life. It may affect accuracy of the results.

24. DISPOSAL AND DESTRUCTION

The test strips are a single-use device and are discarded after each use during the analysis. Since these devices contact with blood, they are infected and/or potentially infected waste (class B). When applied in medical institutions, used devices are discarded in accordance with the sanitary-epidemiological requirements for organizations carrying out medical activities and for this type of waste.

Unutilized devices are referred to class A waste, they do not contain harmful substances and can be discarded with household waste in accordance with SanPiN 2.1.7.2790-10.

In everyday use in domestic conditions used test strips can be disposed of with household waste, after placing them in paper or plastic packaging.

25. WARRANTY OBLIGATIONS

The warranty period is valid for the period of the guaranteed shelf life.

2020 pei-yuan-ming-jen-hsuan-tzu

案號: 000089	號日期: JAN 31 2020
Case No.: 泰博科技股	Date: JAN 31 2020
本文件 份有限公司	之簽名或蓋章, 在臺灣
臺北地方法院所屬民間公證人	北大聯合事務所認證。
	公證人: 徐 暄 詠
Attested at the Bei-Da Partners Notary Public Office	
of Taiwan Taipei District Court, R.O.C., that the	
signature(s)/seal(s) of <u>TaiDoc</u>	
in this document is/are authentic. Technology Corp..	
Notary Public	<u>Rean</u>
	Hsu, Hsuan-Yung

10F-10, No.148, Sec. 4, Zhongxiao E. Rd.,
Da'an Dist, Taipei City 106, Taiwan (R.O.C.)

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