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
Cellprep

The next generation of Liquid Based Cytology

Cellprep AUTO
LBC System
user manual

www.biodyne.asia

 **Biodyne**

Cellprep

Table of Contents

Part.1 Introduction

1-1 The outline of Cellprep	2
1-2 Advantages and Features of Cellprep AUTO System	3

Part.2 Cellprep AUTO LBC System Specification

2-1 Specification and Structure of Cellprep AUTO	5
2-2 Term of each parts of Cellprep AUTO	6
2-3 Touch Screen	8

Part.3 Operating instructions

3-1 Consumables	12
3-2 Check points before operating AUTO	13
3-3 Guide for AUTO Turbidity Smear Mode	14
3-4 Guide for Manual Turbidity mode	22

Part.4 Troubleshooting

4-1 Error Message	30
4-2 How to release air pressure during error	33
4-3 How to treat a vial and a top cap during error	34

Part.5 Maintenance

5-1 Barcode attachment	40
5-2 Cap opener & Suction	41
5-3 slide rack / Alcohol container for cell fixation	45
5-4 Waste Filter Bin	46
5-5 Slide Holder	47
5-6 Waste drainage container	48

Part.1 Introduction

1-1 The Outline of Cellprep

The outline of Cellprep

Cellprep LBC System is a LBC device that smears cells onto slide to make convenient for examination and diagnosis as collecting abnormal exfoliated cells from the human body (from body fluids like urine, ascites etc.). It is compose of Main Frame, Vial to preserve the specimen, Filter, Slide, Brush etc.

Cellprep is developed by Biodyne with the world's first blowing method. Currently, it is technologically developing and manufacturing based on 6 patents within and out of Korea.

Explanation of Cellprep Patents


Liquid-based Cytology (LBC) System for various cancer screenings


Preservative solution for cytological diagnosis


Vial for liquid-based cell preservation and examination

Cervical cancer brush for collecting specimen from vagina

Device that uses air pressure to attach cells

Device that uses blowing method to attach cells

We are continuously on development and research on technology as not only satisfying our current status and also trying and striving to keep pace with the rapid changes in the market.

Cellprep LBC System

Cellprep LBC System is a LBC device and specimens are collected in a Cellprep vial containing Cellprep exclusive solution which helps preserve exfoliated cells from body fluids such as urine, ascetics etc. From these specimens, it smears the cells onto a slide making it easier for diagnosis.

A full kit of Cellprep, to collect mono-layered cell smears for better examination, consists of Cellprep LBC Processor, Cellprep solution vial with solutions for cell preservation, Cellprep membrane filter, Cellprep coating slide and brush.

Through the process, Cellprep LBC System smears the cells on a slide producing a very thin, mono-layered 20mm (diameter) circle of cells. It minimizes misdiagnosis by reducing the amount of overlapping cells on a slide.

Cellprep Solution, exclusively developed by Biodyne, is based on ethanol. It removes unnecessary materials like mucus, inflammatory cells and blood to make examination more convenient. Furthermore, unlike other competing products, Cellprep Solution does not require any pre-preparation. (Exceptionally, additional solution might be necessary in 5~10% of cases depending on the specimen.)


Part.1 Introduction

1-2 Advantages and Features of Cellprep AUTO LBC System

Introduction of Cellprep Auto

Cellprep Auto is a processor that determines the status of specimen automatically and smears as its optimal condition. The full-automatic equipment, Slide and solutions in use at Cellprep 4.7 manually fed automatically through a massive supply one comparison, to match supply.

Cellprep Auto consists of Filter supplier, Cap Opener, Suction, Blower, Slide supplier and solution provider. Filters, vials, vial cap open from the slide supply, the turbidity sensing to measure the cell volume in vials, the cells in vials for implantation in a suction filter, the blower to the cell smear on the slide, the cells are attached to the fixing slide to alcohol the whole process of transferring up to occur automatically, All courses take place making time is very quick.


It is not necessary for user to stay with the device as 40 tests can be processed automatically at once.
Thus, it is effective.

Part.1 Introduction


1-2 Advantages and Features of Cellprep AUTO LBC System

Features of Cellprep AUTO LBC System


Automatic System

Batch processing from solution supply to cell smearing onto the slide.


Multi-Barcode System

QR Code barcode reader matches the patient's sample with the slides and vials, preventing mix-up in advance.


Auto Vortexing System

Auto vortexing disentangles the specimen to improve homogeneity of cells.


Auto Slide Supply & Collection System

Before running, 40 slides are loaded and are collected in a slide tray after smearing.


Auto Capping System

With auto capping system, there's no need to open the cap each time to allow suctioning.


Auto Filter Supply System

50 Filters can be loaded at once and they are supplied automatically.

Auto Vial Supply & Collection System

Before running, 40 vials are loaded and are collected in a vial tray after use.

LCD Touch Screen System

For user convenience pressure and minute settings can be adjusted according to cell concentration and existence of mucus etc.


Solution Specialized for Each Specimen

A certain amount of blood in the specimen collected is naturally dissolved by the solution. (exception of 5-10% of the samples are noted)

Part.2 Cellprep AUTO LBC System Specification

2-1 Specification and Structure of Cellprep Auto

Specification of Cellprep AUTO

Main Features


Rated Voltage / Frequency	AC 220/230V, 50/60Hz
Power consumption	1 KW
Size	1,000(W) x 800(H) x 650(D) mm
Weight	175kg
Operation Environment	Temperature : 15~32°C / Humidity : Under 60% RH
Storage Environment	Temperature : -10 ~ 50°C / Humidity : Under 60% RH


* We are continuously on development and research on technology as not only satisfying our current status and also trying striving to keep pace with rapid changes in the market.

2-2 Term of each parts of Cellprep AUTO

Term of each parts of Cellprep AUTO

No.	Name	Function
1	LCD Control Panel	Controls equipment / shows status
2	Inspection Cover	To check the error of Device
3	Alcohol Rack (Alcohol container for cell fixation)	Load the smeared slide
4	Slide Rack	Load the slide and supply
5	Waste Filter Bin	Collecting the used Filter


No.	Name	Function
6	Filter Magazine	Load the Filter and supply automatically
7	Cap Opener	Fix the vial/ opens and closes the vial cap
8	Suction/Blower Unit	Smear the cells in solution to the slide using air pressure
9	Slide Transfer	Supplying a slide and discharging the smeared slide
10	Slide Holder Unit	Hold the slide for smearing
11	Vial Tray	Load the vial
12	Air combination	Blocking the compressed air from air tank at a predetermined pressure

2-2 Term of each parts of Cellprep AUTO


Term of each parts of Cellprep AUTO

No.	Name	Function
13	Vial Transfer 1	Transferring the vial on vial tray
14	Vial Transfer 2	Transferring the vial secondary
15	Vial Transfer 3	Transferring the solution to the opener
16	Vial Scanning Unit	Barcode recognition unit
17	Vortex	Automatically mix the specimen in the vial
18	Turbidity Defect Sensor	Recognize the density of specimen in the vial and determines the
19	Air Compressor	Compress the air
20	Air Tank	Storing compressed air from the compressor

2-3 Touch Screen

Touch Screen

No.	Functions
1	Display mode according to the amount of cells
2	Specimen mode selection
3	Check the operation status through the status or alerts
4	Values of turbidity measuring
5	Begins Auto run
6	Stops operation
7	RESET ① Revoke the error status ② Move to Initial mode
8	Go to SETUP
9	Move to Manual Setting
10	Reset button (Initialize the device as the default standby)
11	Go to Filter Magazine Front
12	Status window (Check the result of Vail work and Slide work)
13	Shows pressure of Blower
14	Shows pressure of Suction
15	Number of worked slide / Work planned quantity when working 1 cycle
16	Displays number of runs



Part.2 Cellprep Auto LBC System Specification

2-3 Touch Screen

Display is presented in English


User Setup Screen

Please click on the SETUP button for 3 seconds or longer in the main screen.


2-3 Touch Screen

User Setup Screen

① Turbidity


ON : Enable auto-turbidity selection function
OFF : Disable auto-turbidity selection function

② Vortex

ON : Enable vortex function
OFF : Disable vortex function

* During use of auto-turbidity selection function,
Vortex ON will be automatic
Turbidity ON → Vortex ON
Turbidity OFF → Vortex ON/OFF Selection Function


③ Buzzer

ON : Enable warning sound
OFF : Disable warning sound

2-3 Touch Screen

User Setup Screen

④ Barcode

Input of Barcode Digit Number (1~30 digits)


* If you have set 0 as a value, the operation proceeds with barcode function at the OFF mode.

Part.3 operating instructions

3-1 Consumables


Consumables


- ▶ **Cellprep vials** are made from plastic material and the solution within is based on ethyl alcohol. According to the characteristics of the specimens, it is divided into 4 types in order to preserve the morphology of the cells effectively and remove blood/mucus. Also, the cap is made of 2 parts; middle cap and top cap. The middle cap has a Nano polyester filter which can remove impurities and helps with an even cell smear. Vials can be stored up to 2 years in room temperature.
- ▶ **Cellprep Membrane Filter** is composed of a plastic frame, micro membrane and sub membrane.
- ▶ **Cellprep Adhesion(coated) Microscope slide** is coated with a high-quality material and is optimized for Cellprep.
- ▶ It is recommended that lab coats and rubber gloves are worn at all times. Prepare a separate slide rack and alcohol container to insert cell smeared slides right after a cycle is completed.

3-2 Check points before operating Cellprep AUTO


Check points before operating Cellprep AUTO


1. Check if filters are fully loaded in filter magazine

2. Check if slide rack has smeared slide

3. Check if smeared slide is in alcohol container


4. Check if waste filter bin is empty

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18


Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode


a) Turn on the Power.

Press the power button on the bottom right side of the machine to turn on the power.

b) Open the Front cover

Front cover is in the form of a removable, the suction unit and cap opener will be shown when the cover opened.


c) Check the Water Filter

Slightly push the tube on front Suction Unit with hands to discharge the liquid repeating twice and start the operation after compressor stops.


Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode


d) Check Initial Screen

You may disregard this step as 'INITIAL screen' is shown after turning on device.

e) Check Initial Screen


After turning on the power, operate next to Process of INITIAL All Start to make sure all the location of the driver of the equipment is set in default position.

If the INITIAL is not been processed, try again pressing INITIAL All Stop button.

Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode


f) Load the filters in Filter Magazine

Load Filters in the filter magazine in the upper left corner. (Up to 56 filters can be loaded at Once.)

Open the both door on filter magazine. Put the filter in a dedicated rack and facing up the micro membrane of the filter. Pushdown from the top of Filter magazine and gently pull the rack to the left then the Filter fit in and rack only exit. It is better to have sufficient quantity before operate the processor.

In case of the Filter Magazine locate at the backside during operating, press STOP on Touch screen and the magazine will move to the front by pressing Filter MGZ Forward thus, filters can be fill easily.


g) Fill the alcohol in the cell fixing alcohol container

The alcohol container for cell fixation is positioned in the bottom right hand side corner. Simply pour in the alcohol when needed. The container is removable thus if the alcohol is contaminated or when necessary, it may be easily washed.

The slide can be damaged if the position of alcohol container changed, be cautious for its position.

(It is suitable to put give an alcohol according to the height of the intermediate level of alcohol container in the state into a slide rack)

h) Cellprep Adhesion(coated) Microscope slide


Load the Cellprep coating slide, attached bar code to the slide rack.

The bar code attached side should be toward to the Bar code Reader and be cautious for the position of slide rack.

The bar code attached side should be toward to the Bar code reader and be cautious for the position of slide rack


Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode


i) Please click on the **SETUP** button for 3 seconds or longer in the main screen.

The screen will switch to User Setup Screen

j) Select **Turbidity ON**.


Using the auto- turbidity selection

If Turbidity AUTO Mode is selected, a user must check that there is no brush in the vial because the brush in the vial may interfere with the auto-turbidity measurement.


Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode

k) Prepare Cellprep vial


Up to 40 Cellprep vial can be loaded at once. The slide, completed a smear, is come out according to the order of Vial's work.

Put the vials from the front so no vacancy remains.

Check the specimen needs Dilute or Lysis before insert vial.


(● Refer part 6 for guideline of Dilute and Lysis)

Push the tray by aligning the groove on the location of the next section when mounting the vial tray.


Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode


l) Select the types of specimens

Set the mode according to the type of specimen by pressing the buttons on the touch screen.

m) Specify the Quantity for operating

Number press menu comes out when press red number on Working Count.


Enter the desired quantity in the workplace.

Part.3 operating instructions


3-3 Guide for Auto Turbidity Smear Mode

n) Press Reset then press Start button.

Please push 'START' button. The processing of filter in, cap open, suction, blowing and slide transfer is automatically operated. If you want to stop the operation, please push 'STOP' button. Push the START button to continue with the operation. If you want to start from the beginning after stopping the operation, press STOP and then HOME to re-set the position.


If the device does not operate with the START button, please check on the main screen if the message reads 'MACHINE Auto Stop'. (If the mentioned message is not displayed press HOME button to re-set the device and try again. If the screen shows a different message, check its details and resolve the issue.)

o) Check the result of Smear and Fix the slide

When all smearing completes, TOTAL VIAL PROCESSING RESULT displays in the touch screen.

If the work completed in normal, the vial displays with Green, otherwise it displays with Black.

Press the arrow on TOTAL VIAL PROCESSING RESULT or below of SLIDE RESULT, it moves to GLSS PRECESSING RESULT and the SLIDE work can be checked.

Take out the Smearing completed Slide that hold in the cell fixing alcohol container, and fix for 30 minutes in separate 90% alcohol then dye.

215

Part.3 operating instructions

3-3 Guide for Auto Turbidity Smear Mode


p) Arrange the waste liquid in the filter container.

The liquid in the filter container shouldn't be more than 80%, it is better to empty after 1 cycle.

Check the filter container not to overflow.

Maintain not to over 80%, and install not to sink the end of waste liquid line.

(※ There should be extra holes on the lid of filter container to circulate the air.)

q) Turn off the Power.


Press the power button to turn off the power at the bottom right side of the processor.

Part.3 operating instructions


3-4 Guide for Manual Turbidity mode

(From a to h has same turbidity and smear auto mode.)


a) Turn on the Power.

Press the power button on the bottom right side of the machine to turn on the power.

b) Open the Front cover

Front cover is in the form of a removable, the suction unit and cap opener will be shown when the cover opened.


c) Check the Water Filter

Slightly push the tube on front Suction Unit with hands to discharge the liquid repeating twice and start the operation after compressor stops.

Part.3 operating instructions


3-4 Guide for Manual Turbidity mode

(From a to h has same turbidity and smear auto mode.)


d) Press INITIAL button.

You may disregard this step as 'INITIAL screen' is shown after turning on device.

e) Check Initial Screen

After turning on the power, operate next to Process of INITIAL All Start to make sure all the location of the driver of the equipment is set in default position.


If the INITIAL is not been processed, try again pressing INITIAL All Stop button.

28

Part.3 operating instructions

3-4 Guide for Manual Turbidity mode

(From a to h has same turbidity and smear auto mode.)


f) Load the filters in Filter Magazine.

Load Filters in the filter magazine in the upper left corner.
(Up to 56 filters can be loaded at once)

Open the both door on filter magazine. Put the filter in a dedicated rack and facing up the micro membrane of the filter. Pushdown from the top of Filter magazine and gently pull the rack to the left then the Filter fit in and rack only exit. It is better to have sufficient quantity before operate the processor.

In case of the Filter Magazine locate at the backside during operating, press STOP on Touch screen and the magazine will move to the front by pressing Filter MGZ. Forward thus, filters can be fill easily.


g) Fill the alcohol in the cell fixing alcohol container

The alcohol container for cell fixation is positioned in the bottom right hand side corner. Simply pour in the alcohol when needed. The container is removable thus if the alcohol is contaminated or when necessary, it may be easily washed.

The slide can be damaged if the position of alcohol container changed, be cautious for its position.

(It is suitable to put give an alcohol according to the height of the intermediate level of alcohol container in the state in to a slide rack.)

h) Cellprep Adhesion(coated) Microscope slide


Load the Cellprep coating slide, attached bar code to the slide rack.

The bar code attached side should be toward to the Bar code Reader and be cautious for the position of slide rack.

The bar code attached side should be toward to the Bar code reader and be cautious for the position of slide rack


Part.3 operating instructions

3-4 Guide for Manual Turbidity mode


i) Please click on the **SETUP** button for 3 seconds or longer in the main screen.

The screen will switch to User Setup Screen

j) Select **Turbidity OFF**

Disable auto-turbidity selection function


k) **Vortex ON/OFF Selection Function**


Part.3 operating instructions

3-4 Guide for Manual Turbidity mode

1) Prepare Cellprep vial


Up to 40 Cellprep vial can be loaded at once. The slide, completed a smear, is come out according to the order of Vial's work.

Put the vials from the front so no vacancy remains.

Check the specimen needs Dilute or Lysis before insert vial.

(● Refer part 6 for guideline of Dilute and Lysis)


Push the tray by aligning the groove on the location of the next section when mounting the vial tray.

Part.3 operating instructions


3-4 Guide for Manual Turbidity mode

m) Select the types of specimens


Select A, B, C on the touch screen according to the density and condition of prepared specimen.

- **Mode A**
 - Used when amount of collected cells in solution vial are normal
 - Standard mode for operation
- **Mode B**
 - Used when too less amount of collected cells in solution vial
 - Pressure of suction increases to collect general amount of cells to membrane filter
- **Mode C**
 - Used when too concentrated amount of collected cells in solution vial
 - Pressure of Suction decreases to collect general amount of cells to membrane filter

n) Select the types of specimens

Set the mode according to the type of specimen by pressing the buttons on the touch screen.


o) Specify the Quantity for operating

Number press menu comes out when press red number on Working Count.

Part.3 operating instructions


3-4 Guide for Manual Turbidity mode

p) Press Reset then press Start button.

Please push 'START' button. The processing of filter in, cap open, suction, blowing and slide transfer is automatically operated. If you want to stop the operation, please push 'STOP' button. Push the START button to continue with the operation. If you want to start from the beginning after stopping the operation, press STOP and then HOME to re-set the position.


If the device does not operate with the START button, please check on the main screen if the message reads 'MACHINE Auto Stop'. (If the mentioned message is not displayed press HOME button to re-set the device and try again. If the screen shows a different message, check its details and resolve the issue.)

q) Check the result of Smear and Fix the slide

When all smearing completes, 'TOTAL VIAL PROCESSING RESULT' displays in the touch screen.

If the work completed in normal, the vial displays with Green, otherwise it displays with Black.


Press the arrow on 'TOTAL VIAL PROCESSING RESULT' or below of 'SLIDE RESULT', it moves to 'Slide Glass Processing Result' and the SLIDE work can be checked.

Take out the Smearing completed Slide that hold in the cell fixing alcohol container, and fix for 30 minutes in separate 95% alcohol then dye.

Part.3 operating instructions

3-4 Guide for Manual Turbidity mode


r) Arrange the waste liquid in the filter container.

The liquid in the filter container shouldn't be more than 80%, it is better to empty after 1 cycle.

Check the filter container not to overflow.

Maintain not to over 80%, and install not to sink the end of waste liquid line.

(※ There should be extra holes on the lid of filter container to circulate the air.)

s) Turn off the power.

Press the power button to turn off the power at the bottom right side of the processor.

Part.4 Troubleshooting

4-1 Error Message

Message Display	Explanation
Equipment AutoRun Start	Equipment AutoRun is in-progress
Equipment AutoRun End	Equipment AutoRun is at a pause for standby
[A-001] Air pressure low	Air pressure is at a lowered state
[A-002] Left door open	The left door is open
[A-003] Right door open	The right door is open
[A-004] No filters in filter magazine	There is no Filter in the Filter Magazine during AutoRun
[A-005] Not reading a vial barcode	Vial Barcode is not readable during AutoRun
[A-006] Not reading a slide barcode	Slide Barcode is not readable during AutoRun
[A-007] Vial cap not opened at Opener	The top cap of Vial is not opened at Opener Operation during AutoRun
[A-008] Slide glass not loaded	Absence of Slide during AutoRun
[A-009] No Vial at Opener	The Vial is dropped at Opener Operation during AutoRun
[A-010] Problem in filter Magazine or pusher	The operation of Filter Magazin of Pusher Cylinder is in problem during AutoRun
[A-011] Vial tray position not correct (Left/Right)	Vial Tray during a production operation conditions are not in place
[A-012] No slide glass in Loading tray before Autorun	Absence of Slide in Loading Tray before an AutoRun task
[A-013] Slide glass not removed from Unloading tray	The Slides compleated at previous work remain Unloading Tray
[A-014] Same barcode is detected in Loading Tray	Same barcode is detected in Loading Tray
[A-015] Vial cap is not clesed at opener	Vial cap is not clesed at opener

Part.4 Troubleshooting

4-1 Error Message

Message Display	Explanation
[A-101] Motor Alarm - Vial Transfer1 Za	Vial Transfer1 Za Motor / Driver / Sensor is in problem
[A-102] Motor Alarm - Vial Transfer1 Xa	Vial Transfer1 Xa Motor / Driver / Sensor is in problem
[A-103] Motor Alarm - Vial Transfer1 Ya	Vial Transfer1 Ya Motor / Driver / Sensor is in problem
[A-104] Motor Alarm - Vial Transfer3 Za	Vial Transfer3 Za Motor / Driver / Sensor is in problem
[A-105] Motor Alarm - Vial Transfer3 Ya	Vial Transfer3 Ya Motor / Driver / Sensor is in problem
[A-106] Motor Alarm - Filter Transfer Xa	Filter Transfer Xa Motor / Driver / Sensor is in problem
[A-107] Motor Alarm - Slide Transfer Za	Slide Transfer Za Motor / Driver / Sensor is in problem
[A-108] Motor Alarm - Slide Transfer Xa	Slide Transfer Xa Motor / Driver / Sensor is in problem
[A-201] Vial in Opener	The presence of a Vial in the Opener during initialization of the Device
[A-202] Vial in Vial Transfer 2	The presence of Vials in the region of Vial Transfer2 during initialization of the Device
[A-203] Vial in Vial Transfer 1	The presence of a Vial in the Vial Transfer1 Gripper during initialization of the Device
[A-204] Slide in Blower	The presence of a Slide in the region of the Blower during initialization of the Device
[A-205] Vial cap in Opener	The presence of a Vial Cap in the Opener during initialization of the device
[A-301] VT1 Loading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer1 Loading
[A-302] VT1 Unloading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer1 Unloading
[A-303] VT1 Loading2 Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer1 Loading2
[A-304] VT2 Barcode Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer2 Barcode


4-1 Error Message

Message Display	Explanation
[A-305] VT2 Loading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer2 Loading
[A-306] VT2 Unloading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer2 Unloading
[A-307] VT3 Loading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer3 Loading
[A-308] VT3 Unloading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer3 Unloading
[A-309] VT4 Loading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer4 Loading
[A-310] VT4 Unloading Timeover Error	Timeover has occurred due to a problem during operation of Vial Transfer4 Unloading
[A-311] Opener Opening Timeover Error	Timeover has occurred due to a problem during operation of Opener Opening
[A-312] Opener Closing Timeover Error	Timeover has occurred due to a problem during operation of Opener Closing
[A-313] Filter Magazine Timeover Error	Timeover has occurred due to a problem during operation of Filter Magazine
[A-314] Filter Pusher Timeover Error	Timeover has occurred due to a problem during operation of Filter Pusher
[A-315] Filter Transfer Timeover Error	Timeover has occurred due to a problem during operation of Filter Transfer
[A-316] Suction Timeover Error	Timeover has occurred due to a problem during operation of Suction
[A-317] Blower Timeover Error	Timeover has occurred due to a problem during operation of Blower
[A-318] Slide Transfer Loading Timeover Error	Timeover has occurred due to a problem during operation of Slide Transfer Loading
[A-319] Slide Transfer Unloading Timeover Error	Timeover has occurred due to a problem during operation of Slide Transfer Unloading
[A-320] Slide Transfer Barcode Timeover Error	Timeover has occurred due to a problem during operation of Slide Transfer Barcode


4-2 How to release air pressure during error

How to release air pressure during error

If you want to release the air pressure during error condition, rotate a red valve above the suction unit

Now you can move each part of the equipment by hand after all air pressure in the machine has been removed


Do not put your hand inside the operating machine

Part.4 Troubleshooting


4-3 How to treat a vial and a top cap during error

(if the device is stopped and if a vial and a top cap are left)

Troubleshooting -1

a) Click the All Gripper Open button in the Initial Screen.

b) After Click, open the both doors with your hands and place all the vials in order.

Place the doors next to both sides of the equipment (see the arrow mark)

If the vial in error or the top cap still remains on the cap opener although you have tried so, you will have it removed according to <Troubleshooting -2 in the case where the equipment is stopped during operation with the vial or the top cap remaining > on the next page.

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
Do not put your hand inside the operating machine

Part.4 Troubleshooting


4-3 How to treat a vial and a top cap during error

(if the device is stopped and if a vial and a top cap are left)


Troubleshooting -2

a) Click the Manual button on the Main screen for 3 seconds or more.

b) Click the User Manual 2 (Opener) button on the User Manual 1 (Cleaning) screen.

c) Click the Opener Ready button on the User Manual 2 (Opener) screen.


Part.4 Troubleshooting


4-3 How to treat a vial and a top cap during error

(if the device is stopped and if a vial and a top cap are left)


Troubleshooting -2

d) Click the Opener Y1 Forward button.

Opener Y1 moves forward

e) Click the Opener Gripper Open button.


Operators must be trained before they are allowed to carry out hazardous work

Part.4 Troubleshooting


4-3 How to treat a vial and a top cap during error

(if the device is stopped and if a vial and a top cap are left)


Troubleshooting -2

The vial hanging in Opener Y1 Gripper is removed
(Take caution because it may fall abruptly.)

1) Click the Opener Y2 Forward button.

Opener Y2 has come out forward.
Remove the vial's primary filter.


Operators must be trained before they are allowed to carry out hazardous work

Part.4 Troubleshooting


4-3 How to treat a vial and a top cap during error

(if the device is stopped and if a vial and a top cap are left)


Troubleshooting -2

g) Click the Opener Y1, 2 Backward button.

Opener Y1 and Y2 all go back to their original positions.

h) Click the Opener End button.


 Operators must be trained before they are allowed to carry out hazardous work

Part.4 Troubleshooting

4-3 How to treat a vial and a top cap during error

(if the device is stopped and if a vial and a top cap are left)

Troubleshooting -2

i) Click the Initial All Start button in the Initial Screen.


Return to the Main screen.

5-1 Barcode attachment


Barcode attachment

One-dimension bar code and QR code can be used, but a Micro QR code is not supported.

Notes

1. The bar code of the vial must be attached within 35mm as shown in the above figure. (Especially, don't forget to have it attached vertically in the case of one-dimension barcode.)
2. Be careful not to have the adhesive side of the barcode revealed to outside.
3. For the QR code, the standard size is horizontal (10 ~ 12mm) x vertical (10 ~ 12mm)
4. The bar code must be printed with a margin of 0.5mm or less vertically and horizontally



Notes

1. The Slide bar code must be attached to a size of 25mm or less, appropriate for a white printed area as shown in the above figure.
2. For the QR code, the standard size is horizontal (7~8mm) x vertical (7~8mm)
3. The bar code must be printed with a margin of 0.5mm or less vertically and horizontally

5-2 Cap opener & Suction


Cap opener & Suction

a) Click the Manual button on the Main screen for 3 seconds or more.

5-2 Cap opener & Suction


Cap opener & Suction

b) ①Click the **Cleaning Ready** button.

Initialization is in progress


After initialization, Cleaning Ready will be changed to Cleaning Start while blinking

c) ②Click the **Opener Forward** button.

Opener Y1 and Y2 Cylinder have come forward.

Wipe the Opener with cloth soaked with alcohol, and pour about 20 ml of alcohol (70-95%) into the opener.

d) ③Click the **Opener Discharge** button.


Confirm that the alcohol poured into Opener has disappeared.

Operators must be trained before they are allowed to carry out hazardous work


5-2 Cap opener & Suction

Cap opener & Suction


e) ④Click the Opener Backward button.

Opener Y1 and Y2 Cylinder have come back to their places

f) ⑤Click the Suction Up button.

Using the wash bottle and laboratory instruments, etc, Pour about 20 ml of alcohol (70-95%) into the Suction Unit

g) ⑥Click the Suction Discharge button.


Confirm that the alcohol poured into the suction unit has disappeared.

Operators must be trained before they are allowed to carry out hazardous work


5-2 Cap opener & Suction

Cap opener & Suction


h) ⑦ Click the Suction Down button.

The suction unit is being lowered


i) ⑧ Click the Cleaning End button.

j) Click the Initial All Start button in the lower left of the Initial Screen to return to the Main screen.

5-3 Slide rack / Alcohol container for cell fixation

Slide rack / Alcohol container for cell fixation

① Residues of the glass-dust may be visible on the slide rack. Use gauge to wipe after operating the device.

② For the alcohol rack to fix cells, have 95% alcohol filled up to the line before operating equipment.


95% alcohol is recommended to be replaced about one or two times a day in order to avoid contamination

After completing a day's jobs, remove all of the alcohol contained in the alcohol rack for fixing cells, wash it with flowing water, and gently wipe it clean with gauze

5-4 Waste Filter Bin

Waste Filter Bin


The liquid in the filter container shouldn't be more than 80%, it is better to empty after 1 cycle.

Check the filter container not to overflow. Maintain not to over 80%, and install not to sink the end of waste liquid line.

5-5 Slide-Fixing Unit (Slide Holder)

Slide-Fixing Unit (Slide Holder)


After completing a day's jobs, gently wipe with soft cloth to prevent glass fragments, liquid, and others from being left on the slide.

52

5-6 Waste drainage container

Waste drainage container

Please check the level of waste liquid before and after operating the device. Maintain the level less than 80% of the container.

In case of that the location of waste drainage bottle is too high or the drainage hose is twisted or the bottle is filled with liquid waste, the liquid waste cannot be discharged well.

User can use any kinds of bottles for drainage, not supplied bottle. But the pathway (hose) in the bottle is required for air to escape.

Cellprep

Biodyne Co., LTD

#301, Kranz Techno Bld., Dunchon-Daero 388 beon-Gil, Jungwon-Gu, Seongnam-si, Gyeonggi-Do, Korea (13403)

Tel) +82-31-741-7253 / Fax) +82+31-741-7259

E-mail) info@biodyne.asia

www.biodyne.asia

For safe use of device, please be sure to read instruction correctly.

Pictures on this user's guide may differ to the actual shape of the product part

Appendix №1 to the User Manual
for Cellprep AUTO LBC System
for Russian Federation only

The User Manual for **Cellprep AUTO LBC System** developed by **Biodyne Co., Ltd., Korea** in accordance with inner and international requirements is standard for all countries. In compliance with local requirements and to observe the Russian legislation, the User Manual supplied along with the instrument to Russian market will be supplemented with following additional data. Data contained in this Appendix shall be considered as prioritized.

1. The following information is presented in current and new Subsections into "Part.1 Introduction. 1-1. The Outline of Cellprep" as shown below


Overview of Cellprep

Cellprep AUTO System is a LBC (Liquid Based Cytology) device. Cellprep AUTO System is an auxiliary equipment intended for preparation of a mono-layered cytological sample using LBC method for its further staining and examination with the microscope for the purpose of early diagnosis of precancerous conditions, detection of oncopathology and the degree of its malignance.

Specimens are collected in special Cellprep vials containing exclusive solution which helps preserve exfoliated cells from body fluids such as pleural effusion, pericardial effusion, ascitic fluid, sputum, urine, cerebrospinal fluid, synovial fluid and Broncho alveolar lavage (bronchial washings). From these specimens, it smears the cells onto a slide making it easier for diagnosis.

Cellprep is developed by Biodyne with the world's first blowing method. Currently, it is technologically developing and manufacturing based on 6 patents within and out of Korea.

Explanation of Cellprep Patents

We are continuously on development and research on technology as not only satisfying our current status and also trying and striving to keep pace with the rapid changes in the market.

Cellprep LBC System

A full kit of Cellprep, to collect mono-layered cell smears for better examination, consists of Cellprep LBC Processor, Cellprep vial with solutions for cell preservation, Cellprep membrane filter, Cellprep coating slide and brush.

58

**Appendix №1 to the User Manual
for Cellprep AUTO LBC System
for Russian Federation only**

The User Manual for **Cellprep AUTO LBC System** developed by **Biodyne Co., Ltd., Korea** in accordance with inner and international requirements is standard for all countries. In compliance with local requirements and to observe the Russian legislation, the User Manual supplied along with the instrument to Russian market will be supplemented with following additional data. Data contained in this Appendix shall be considered as prioritized.

1. The following information is presented in current and new Subsections into "Part.1 Introduction. 1-1. The Outline of Cellprep" as shown below

Overview of Cellprep

Cellprep AUTO System is a LBC (Liquid Based Cytology) device. Cellprep AUTO System is an auxiliary equipment intended for preparation of a mono-layered cytological sample using LBC method for its further staining and examination with the microscope for the purpose of early diagnosis of precancerous conditions, detection of oncopathology and the degree of its malignance.

Specimens are collected in special Cellprep vials containing exclusive solution which helps preserve exfoliated cells from body fluids such as pleural effusion, pericardial effusion, ascitic fluid, sputum, urine, cerebrospinal fluid, synovial fluid and Broncho alveolar lavage (bronchial washings). From these specimens, it smears the cells onto a slide making it easier for diagnosis.

Cellprep is developed by Biodyne with the world's first blowing method. Currently, it is technologically developing and manufacturing based on 6 patents within and out of Korea.

Explanation of Cellprep Patents

Liquid-based
Cytology (LBC)
System for various
cancer screenings

Preservative solution
for cytological
diagnosis

Vial for liquid-based
cell preservation
and examination

Cervical cancer brush
for collecting specimen
from vagina

Device that uses air
pressure to attach
cells

Device that uses
blowing method to
attach cells

We are continuously on development and research on technology as not only satisfying our current status and also trying and striving to keep pace with the rapid changes in the market.

Cellprep LBC System

A full kit of Cellprep, to collect mono-layered cell smears for better examination, consists of Cellprep LBC Processor, Cellprep vial with solutions for cell preservation, Cellprep membrane filter, Cellprep coating slide and brush.

Through the process, Cellprep AUTO LBC System smears the cells on a slide producing a very thin, mono-layered 20mm (diameter) circle of cells. It minimizes misdiagnosis by reducing the amount of overlapping cells on a slide.

Cellprep Solution, exclusively developed by Biodyne, is based on ethanol. It removes unnecessary materials like mucus, inflammatory cells and blood to make examination more convenient. Furthermore, unlike other competing products, Cellprep Solution does not require any pre-preparation. (Exceptionally, additional solution might be necessary in 5~10% of cases depending on the specimen.)

Authorized representative responsible for technical service and support in Russia

BioLine LLC
197046, Saint-Petersburg,
Pinsky alleyway, 3, lit. A
Phone: +7(812) 320 49 49
Hotline: +8 800 333 00 49
Web: <http://www.bioline.ru>

Indications and contraindications to the use, adverse events

Indications for use: preparation of a mono-layered cytological sample using LBC method for further cytological examination for the purpose of early diagnosis of precancerous conditions, detection of oncopathology and the degree of its malignance

Contraindications: no

Adverse events: no applicable.

Information about potential consumers of a medical product

Only for the professional users.

Cellprep AUTO is intended for use only by qualified laboratory personnel (for example, a doctor of clinical laboratory diagnostics or medical laboratory technician). All authorized laboratory personnel to work with this device must read the relevant User Manual and should be familiar with the technical features of the instrument before switching for the first time.

The type of test samples

The test samples by Cellprep AUTO: biopsies, biological tissues and body fluids, namely: pleural effusion, pericardial effusion, ascitic fluid, sputum, urine, cerebrospinal fluid, synovial fluid, Broncho alveolar lavage (bronchial washings), except of the blood.

2. New subsections are added into the Part.2 "Cellprep AUTO LBC System Specification."

2-1.Specification and Structure of Cellprep AUTO" and presented as follows

	Main Features	
	Rated Voltage/Frequency	220-230 V ~, 50/60Hz

General Spec.	Power consumption		Max 1 KW
	Size (W*D*H)		1,000 mm X 650 mm X 800 mm
	Overall size allowance		± 5 mm
	Weight (kg)		175 Kg
	Weight allowance		± 5 kg
	Classification by Standard IEC 61010-1:		Class of protection 1 Pollution degree 2 Overvoltage category II Class equipment A
	IP according to IEC 60259		IPX0 (standard)
	Noise measurement		Less 64.5 dBA
	Protection against excess power		7 circuit protectors: - CP4 circuit protector SERVO POWER - CP3 circuit protector DC POWER - CP2 circuit protector PLC POWER - CP1 circuit protector SMPS POWER - MCCB1 circuit protector Power of transformer and compressor - M/C circuit protector Magnetic Contactor - ELB circuit protector Main POWER See more information in section 3-2-2.
	Material and the corrosion resistance		- Cold Steel Plate Steel(SPC) - Galvanizing (0.04~0.08) mm & Powder coating (0.1~0.2) mm - white/grey
	Operation Environment	Temperature	15~32 °C, special manufacturer requirements
		Humidity	20~60% RH well-ventilated premises
	Storage and transport Environment	Temperature	-10~50 °C, special manufacturer requirements
		Humidity	20~60%RH
	Operating altitude		not more 2000 m above sea
	Filter		Max. 56 ea

	Storage Capacity	Slide	Max. 40 ea
		Vial	Max. 40 ea
Instrument	Operation Method		Air compressor, Electric motor
	Compressor		Miniature Compressor (2 ea, 60 LPM)
	Filter Supply Method		Magazine (2 X 28EA)
	Cap Open	DC reducer motor	
		Timing Pulley & belt	
		3 Pin rotation Top Cap Open & Closed	
	Draining		Water Pump (Opener, Suction, Suction Filter, 2.5LPM)
	Suction		Ejector for Suction (12 LPM)
	Slide Holder		Using Slide Transfer – Holding to prevent contamination
	Vortex		<ul style="list-style-type: none"> Allowable weight to shake: 100 g
Electric	Multi-Barcode System		<ul style="list-style-type: none"> LED scanner, excluded from the group of risk (LED safety in accordance with IEC 62471:2006) High-speed 2D scanning Editing function LED targeting Direct light source in the scanner: 1 green LED Scanning light source: 2 red LEDs
	Software		Version XG5000 Ver. 4.06 is specified at the top right of the main screen
	The security class of software in accordance with the Russian standard ГОСТ Р МЭК 62304-2013		Class A. Since the installed software is only the controlling software / firmware of the equipment, no injury or damage to health is possible. No crushing hazard could occur. Also the firmware does not make any clinical decisions.
	Setting and Control		PLC
	Display and Interface		7" TFT Color Touch Screen
Pressure	Pressure of water pump		Max Pressure: 160kPa Max Vacuum: 70kPa

	Pressure of air compressor	Max Pressure: 700kPa Max Vacuum: 68kPa
	Pressure of air tank	Max Pressure: 785kPa
Modes of operation	<u>The setting time of the operating mode and duration</u> - for automatic mode - for manual mode	- during the installation of the System Cellprep AUTO by the operator, the setting time of the operating mode is no more than 60 minutes. Daily setting of the automatic mode by the operator is not required. - during daily operation, the operating mode is set after 1 min. after switching on. Duration of one run is (45±10) min. The operating mode of the System – continuous, not less 8 hours. - not applicable

Storage and transport conditions of Cellprep AUTO

1. Temperature and humidity: -10 °C to + 50 °C special manufacturer requirements; (20-60) % RH, no condensation
2. It must be stored in a place it will not be splashed with water.
3. Barometric pressure, temperature, humidity, ventilation, sunlight, dust, salt, to be stored in a fear-free place of occurrence of adverse effect by such an atmosphere containing such as sulfur content.
4. Tilt, vibration, shock (including during transport), and the like to note stable.
5. Do not store to the occurrence location to the storage location and gas chemicals.

3. The information of Cellprep AUTO standard delivery is added into the Part.2 "Cellprep AUTO LBC System Specification. 2-2.Term of each parts of Cellprep AUTO" and presented in new subsection as shown below

Cellprep AUTO LBC System standard delivery

Cellprep AUTO LBC System standard delivery includes:

1. Cellprep AUTO LBC System – 1 pcs.
2. Alcohol Rack (Alcohol container for cell fixation) – 1 pcs.
3. Slide Rack – 1 pcs.
4. Waste Filter Bin – 1 pcs.


5. Vial Tray – 2 pcs.
6. Main power cable – 1 pcs.
7. User Manual – 1 pcs.

4. The title of section “Part 3: Operating instructions. 3-1 Consumables” is changed and presented as shown below:

Part 3. Operating instructions

3-1. Devices and materials used with Cellprep AUTO

Consumables

- These devices and materials used with Cellprep AUTO are not included in the delivery kit and are purchased by the user on its own.
- **Cellprep vials** (certificate №ФС3 2010/07117 points 5, 6, 7, 8) are made from plastic material and the solution within is based on ethyl alcohol. According to the characteristics of the specimens, it is divided into 4 types in order to preserve the morphology of the cells effectively and remove blood/mucus. Also, the cap is made of 2 parts; middle cap and top cap. The middle cap has a Nano polyester filter which can remove impurities and helps with an even cell smear. Vials can be stored up to 2 years in room temperature.
- **Cellprep Membrane Filter** (certificate №ФС3 2010/07117 p. 10) is composed of a plastic frame, micro membrane and sub membrane.
- **Cellprep Adhesion (coated) Microscope slide** (certificate №ФС3 2010/07117 p. 11) is coated with a high-quality material and is optimized for Cellprep.

- Contact the authorized representative in your country for availability of the consumables.
- It is recommended that lab coats and rubber gloves are worn at all times. Prepare a separate slide rack and alcohol container to insert cell smeared slides right after a cycle is completed.

5. Information on placement, installation of the equipment, instructions for connecting to power supply and protective ground is added to the section: "Part 3. Operating instructions. 3-2. Checks points before operating Cellprep AUTO" and is presented as shown below:


3-2-1. Installation site

1) The equipment should be placed at distance of at least 20 cm from the wall or other device and at least 50 cm from the heating devices, to ensure free ventilation of air around the instrument and to prevent its overheating.

Min. size of the installation site 1200 x 1000 x 850 mm.

2) To ensure a stable operation, the instrument should be placed on an even horizontal surface. The floor of the working area should be even, the table on which the instrument is placed must have a reinforced structure and be strong to withstand a load of 200 kg / m².

3) It is also necessary to provide a place for a waste drainage container. The container is usually placed under the table on which the instrument is.


4) Installation and setup of the equipment should be carried out only by a trained engineer.

3-2-2. Power Supply

Use the specific power cable to connect the equipment to power supply. If another cable is used, electric energy input may change, which will lead to unstable operation of the equipment.

Delivery kit includes a power cable with Schuko plug.

For Cellprep AUTO, it is recommended to use a grounded outlet. DO NOT connect the instrument without ground connection.

Before inserting the plug into a power outlet in the laboratory, make sure that the power supply voltage of the device is suitable for the power line of the laboratory. Connecting the device to the wrong power sources will cause its damage.

The power outlet to which the instrument is connected must be visible and easily accessible to the operator, and should only be used to connect the Cellprep AUTO.

To protect against overcurrent, the instrument is equipped with circuit breakers:


No.	Image	Description
1		ELB (32KGRd/15A) ELB (Earth Leakage Breaker) prevents leakage of power, which can cause damage to electric cords or moisture intrusion. When the switch is up, the Trip Coil Turns on the Breaker. Trip Coil: Coil cable that turns on if the Breaker is turned On. It operated other components by supplying electricity.
2		M/C (MC-12b/AC110V) M/C (Magnetic Contactor) controls electric currents. Large Electric currents are connected or blocked, and smaller Electric currents are controlled.
3		MCCB (BKM-22P/20A) MCCB (Mold Case Current Breaker) protects the machine circuit. It prevents over currents and short circuit currents.
4		CP1, CP2, CP3, CP4 (GCP-32AM-20A/10A/5A) CP (Circuit Protector) helps the Breaker. It protects the circuit by playing the role of a Fuse.



3-2-3. Preparation for starting the instrument

Before starting the work, while the instrument is in a safe condition, it is necessary to verify visually that the instrument is placing evenly on a horizontal surface, as well as to check the instrument and the power cable for external damages.

- Connect the power cable to the socket - ① on the rear left part of the equipment main body
- Connect the waste liquid tube to the hole - ② on the rear left part of the equipment main body
- After the installation is completed, place a waste drainage container under the equipment.
- Make sure that the tube does not reach the bottom of the waste drainage container.
- Insert the plug into a power outlet.
- Move the power button - ③ to I position (ON)



		
<p>① Main power cable</p>	<p>② Hole for waste liquid tube connection</p>	<p>③ Power button</p> <p>I position – ON</p> <p>O position - OFF</p>


- Adjust the settings and confirm the numbers
- Check operating sound of the compressor. Only begin operation when the sound is mute.
- Load in membrane filters by opening up the filter magazine
- Load in the vial containing the specimen
- Insert the slide into the slide tag
- Check the following points before operating


1. Check if filters are fully loaded in filter magazine

2. Check if slide rack has smeared slide

3. Check if smeared slide is in alcohol container.

4. Check if waste filter bin is empty

- Push the START button

If the equipment is operated in violation of the manufacturer's requirements, its protection level may be reduced.

6. Subsection "3-5. Safe instrument switching-off" is added to section "Part 3. Operating instruction".

- 1) To switch off the instrument safely, move the power button on the lateral side of the instrument to the O position ("OFF").
- 2) The system is switched off.
- 3) Disconnect the instrument from the power supply.
- 4) Disconnect the accessories, if necessary.

7. Clearing recommendation, warnings, the service information, responsibility of the Manufacturer and warranty are included into the Part.5 "Maintenance" and shown below

During service, cleaning and disinfection of the instrument, it is recommended to use always laboratory coats and protective gloves.

5-7. Cleaning and storage

- 1) Switch off the power supply.
- 2) Plug off the device.
- 3) It is recommended to use alcohol (not less than 63%) disinfecting agents (e.g., SEPTOCIDE R PLUS, AKHDEZ RED), ethanol 70% or any other disinfectants with corresponding characteristics.
- 4) When cleaning the external surfaces of the instrument, it is not recommended to take off the side panels.
- 5) Do not allow disinfectants or detergents to contact with the scanner or sockets.
- 6) Wipe or spray disinfectant on the external surfaces of the instrument and leave it for 5-10 minutes
- 7) Use a suitable detergent (soapy water or liquid detergent) to clean the surfaces.
- 8) Dry with paper or a clean cloth.
- 9) Regularly check the liquid level in the waste drainage container and, if necessary, empty it.
- 10) Drain out the hose by inserting alcohol into the Suction Unit to prevent clogging, see section 5-2.
- 11) For servicing and replacing the hoses connected to the cylinders, you should contact the service department, see section 5-10.
- 12) Perform regular checkups, cleaning and repair.
- 13) Recommended storage and method.
 - Place where it is water-free.
 - Place where it is dust and humidity free.
 - Place where it is free from slope, vibration, impact.
 - Place where there's no direct sunlight.

5-8. Warnings

Warning during operation

- 1) Allow to work only trained personnel of User Manual and an instruction on labor protection
- 2) Do not reuse solution vial, filter, slide and brush.
- 3) Fixate the solution vial, filter and slide properly.
- 4) After-use slide is automatically transferred to the waste filter bin. Check if there's enough solution to be immersed in. (Enough amount to immerse smeared cell in size of 20mm)
- 5) All laboratory facilities should be complied with local fire safety requirements and have the fire -fighting equipment.

- 6) Characteristics of flammable substances and materials used in the work:

Ethanol is a colorless, flammable liquid. The boiling point +78.3 ° C, flash point + 13 ° C, the auto ignition temperature + 400 ° C. Inflammation limits: lower + 11 ° C, upper + 41 ° C.

- 7) To reduce the ethanol handling risk it is recommended the alcohol container for cell fixation not to fill above the indicated line. After completing a day's jobs, remove all of the alcohol contained in the alcohol rack for fixing cells, wash it with flowing water, and gently wipe it clean with gauze.

Daily norm of alcohol needs in Cellprep AUTO is less 2,0 liters. Ethanol used for work should be stored in accordance with local requirements.

For Russian Federation in accordance with the rules implemented by ПНД Ф 12.13.1-03 «Guidelines. Safety in working in analytical laboratories (general provisions) ».

- 8) In case of fire, immediately disconnect the Cellprep AUTO from the power source. Recommend to extinguish with fine water or fire extinguishers.
- 9) Violation of the operating rules of equipment provided by the manufacturer may be the reason for the deterioration of the protection applied in this system.

Warning during breakdown and remodeling

- 1) When the device is malfunctioning, certified technician must be in charge of the repair.
- 2) Perform regular checkups and repairs in line with monthly inspection.
- 3) When the machine was not-in-use state for long time, perform operation after inspection.
- 4) Do not disassemble or remodel without supervision.
- 5) Stop the operation when malfunctioning and contact the Service Center at Biodyne Co. Ltd. Tel: +82-31-741-7253 or contact to authorized representative in Russia **BioLine LLC**, tel: +7 (812) 320-49-49; 8-800-333-00-49.
- 6) Inspection of equipment after repair is carried out by the Service Center of the company Biodyne Co., Ltd. or an authorized representative.

Warnings, precautions and / or measures taken in the case of external factors such as external electromagnetic fields, electrostatic discharges, radiation, atmospheric pressure and its differences, humidity and air temperature, which can affect the functioning of a medical device.

- 1) Recommended to place the system in places where the intensity of the external electromagnetic radiation and the electrostatic discharge is below the maximum permissible intensity.
- 2) Do not use this equipment near sources of the strong electromagnetic radiation, because it may interfere to its normal functioning.
- 3) Recommended to use the device in the specified operation environmental conditions (see Section 2).

Warnings, precautions and / or measures taken in the case of a predictable risk of electromagnetic interference from a medical device to other medical devices and equipment.

The level of electromagnetic interference of Cellprep AUTO does not exceed the permissible and not have a significant impact on the operation of the other medical devices or an equipment.

5-9 Responsibility of the Manufacturer

Biodyne Co., Ltd., Korea is the legal manufacturer that is responsible for the quality, efficiency and safety of the registered medical product for use, normal operation and possible negative results of using them, except for cases of negligence, or incorrect use through the fault of user, or force majeure.

5-10. The repair and service information

Maintenance and current repairs of devices manufactured by Biodyne Co., Ltd. is available anywhere in the world in one of the international service centers, as well as in your area. To arrange service or repair of your instrument, please contact your nearest office.

For service organizations and technical support see details below.

Service Center

Tel.: +82-31-741-7253

E-mail: info@biodyne.asia

For the Russian Federation the authorized representative:

BioLine LLC

Pinsky lane 3 letter A,
197046, Saint Petersburg,
the Russian Federation

E-mail: service@bioline.ru

Tel: (812) 320-49-49

Fax: (812) 320-49-40

Service Hotline: 8-800-333-00-49.

This instrument is an in vitro medical device and it is usually installed in a laboratory where specimens and other biological tissues are present. For your safety, it is therefore required to clean and disinfect before entering in contact with it, refer to section 5-7.

It is required in case of return that the unit is cleaned and disinfected before sending it back to Biodyne Co., Ltd.

International health rules require that shipments of biohazard materials are do not in standard packages (risk of sanctions).

5-11. Warranty and the service life of device

Warranty applies to medical instruments manufactured by Biodyne Co., Ltd.

Biodyne Co., Ltd. guarantees that the contractual product delivered has been subjected to a comprehensive quality control procedure based on international standards, and that the product faultless and complies with all technical specifications.

The warranty period of storage is not less than 6 months.

The warranty period of operation is 12 months from the date of shipment.

This warranty becomes invalid in case of mishandling, damage and negligence and further in case of usage of inappropriate spare parts and/or accessories or unauthorized modification of the unit.

The instrument service life is 7 years if you implement appropriate maintenance and inspection (self-certification).

Using of equipment in other than specified ways may void the warranty provided by the manufacturer.

8. Information regarding Cellprep Solutions, Gynecology sample preparation, Non gyn sample preparation, Recommended Staining Method of Cellprep and A. B. C. Suction Numerical standard values are presented in new corresponding sections and shown as follows


Part. 6. Cellprep Solution

6.1. Introduction of Cellprep Solutions

Cellprep solution is practicable to preserve specimen in vial for two years as it works for the pre-fixation and preservative of collected samples and is specified on the label. The specimen contained within the vial can be tested up to 3 months of collection.

After processing of pre-preparation with collected cell manually, user can put them into Cellprep vial solution

- The collected specimen should be left in the vial for at least 30 minutes for fixation before running in the Cellprep device for smear.
- But, in case it is difficult to fix the smeared sample in 95% Ethanol instantly, it is possible to use Cellprep Plus Solution (pre-fixation) in Cellprep for more than 30 minutes.
- In Case of some specimens such as Body Fluid, Bronchial Washing, CSF, it has to preserved in Cellprep Vials for 12hours in order to have good quality of cell morphology.


Specialties of Solution

- Great at removing red blood cell, mucus, and impurities without cell transformation
- Possible to run Generic examination

Above Cellprep Solutions are appropriate for its specialties as considering its specimen and possible to preserve cells for long period of time without cell transformation. It is exclusively for Cellprep LBC Processor presenting great staining for better and convenient examination.

Vials of Cellprep Solutions are developed as right for Cellprep LBC Processor to use world's first blowing technology to smear cells on slide.

Cellprep has two kinds of filters to reduce misdiagnosis during the examination caused by impurities. Attaching Cap Primary Filter to vial to make possible to filter impurities within collected specimen, which irritates examination for the first time then secondarily use Membrane Filter to smear pure cells on slide.

Filters on Middle Cap

6.2. Introduction of Solutions for Pre-preparation

Cellprep Solution for Pre-preparation (Not a compulsory process)

Solutions in Cellprep vial remove blood, mucus and some impurities. This is not a compulsory process but in order to have good quality, pre-preparation is recommended before performing the test for specimens that contains a lot of blood appearing dark red or for specimens that looks milky and opaque due to excess mucus.

Pre-preparation Solution		Pre-preparation Solution	Usage
	To dilute the sample. It is commonly used for both Gyn & Non-Gyn (In case of high concentration or too much cells, add about 10~20ml)		To remove irritant elements like mucus and red blood cells. It is divided into two solution; Gyn and Non-Gyn Solution
Dilute Solution		Lysis Solution	

Specimens which requires pre-preparation

- In case of a lot of blood in sample
- In case of a lot of mucus in sample
- In case of sample which has a lot of impurities or high concentration

Introduction of Solutions for Pre-preparation

(1) Lysis

1. Pour a certain amount of specimen into 50ml centrifuge tube and vortex for 10 seconds then centrifuge in 1800rpm for 5min.
2. Remove supernatant and add 20ml of Cellprep Lysis Solution then vortex for 1 minutes and centrifuge in 1800rpm for 5min.
3. Remove supernatant and put sediment into Cellprep Vial.

(2) Dilute

If specimen doesn't smear well due to too much cells or mucus, run operation after adding 20ml of dilute solution and vortexing for 30sec to loosen up.

Part 7. Gynecology sample preparation

7.1. Gynecology System

Cellprep Solution for Pre-preparation (Not a compulsory process)


■ Specimen Collection and precautions

➤ Specimen collection method

1. Adhere the tip of the brush to the endo cervix and rotate the brush in one direction for 2 full turns.
2. Insert the brush in the Cellprep vial and rinse the brush 5~10 times so that the cells are well dispersed.
3. If there seems to be plenty of cells inside the vial, discard the brush. If not enough cells are collected, remove the brush tip and leave it inside the vial.
(※ But, if the brush used is not a cervical brush (as shown in the first picture on the left) always leave the separate brush tip within the vial.)

➤ Precaution

1. Cervical brushes are the most adequate brush to collect specimens. In certain cases, a cyto brush and a plastic spatula can be used together but a cotton swab cannot be used.
2. It is advisable to rotate the brush twice. If less than 2 turns there might be insufficient cells for examination and if it exceeds 2 turns bleeding may occur.
3. Pre-preparation for Gyn: less than 5%, in case of needs

7.2. General Gyn processor

General Gyn processor

1. Put the collected sample into Cellprep Vial
 - Insert the brush in the Cellprep vial then shake it few times until specimen fall off from the brush then remove it from vial.

(Please remove the brush from vial before operate turbidity auto mode with Cellprep Auto processor)

2. Before insert the Cellprep Vial which included specimen to Cellprep LBC processor, please pre-prepare enough for more than 30 minutes.
 - Vortex or shake the Cellprep Vial until specimen fall off when you directly see with naked eyes.
3. Label the slide and prepare with Cellprep LBC Processor.
 - Set up the mode depending of the amount of cells.
(A=Regular, B=Small, C=Overload).
4. After smearing insert the slide into 95% ethanol immediately and let it fix for 30 min. Follow Papanicolaou staining method to stain the slides.


Part. 8. Non gyn sample preparation

8.1. Urinary System (kidney, ureter, urinary bladder, urethra)

- Collecting specimen and pre cautions : Store in a Refrigerator for 2 hours (Before putting urine into Cellprep vial)
- ☐ Voided urine: Take mid-stream urine (Recommended to take urine over half amount of urine cup)
- ☐ Catheterized urine: It is collected by catheter which insert it to urethra. This is for the patient who has difficulty with Voided urine or who has lesions at lower part of urethra.
 1. Let patient drink water every 30 minutes for 3 hours
 2. Intubate until the suspicious lesion part.
 3. Throw away the gathering urine through intubation after 2 hours.
 4. Transfer the urine, which received after 3 hours, to pathological laboratory and centrifuge then put into Cellprep vial.
- ☐ Bladder washing: When cystoscopy (bladder examination), after wash with saline or ringer, take the washing liquid.

■ General Urine processor.

Fig1. Collecting sample	Fig2. Centrifuge	Fig3. Pouring to Cellprep vial, Vortex	Fig4. Insert Cellprep vial to Cellprep LBC Processor
			
	Tube(50ml)/1800rpm, for 5mins	Collect sediment (Cellprep vial)	Shake 4-5 times and insert

1. Pour collected sample into 50ml centrifuge tube and centrifuge in 1800rpm for 5 minutes. (Fig.2).
2. Remove supernatant then pour the sediment (less than 1ml) into Cellprep Vial and leave to fix for more than 30 minutes (Fig. 3).
 - In case of a lot of blood and mucus in sediment, use lysis solution then process next step.
3. Label the slide and prepare with Cellprep LBC Processor(Fig.4).
 - Before insert Cellprep Vial to Cellprep LBC System, Vortex for 5~10 seconds or shake 3~4 times.
4. Fix the slide directly in 95% Ethanol (more than 30 minutes), then stain with Papanicolaou staining for examination.

8.2. Cerebrospinal fluid

- Collecting specimen and precautions: Should not be stored in a Refrigerator for more than 2 hours (before putting cerebrospinal fluid into Cellprep vial.)
 - Cerebrospinal fluid is a clear colorless bodily fluid that occupies the subarachnoid space and the ventricular system around and inside the brain and spinal cord. The usual collection amount is a very small and because it is collected while in lumbar anesthesia so exception care is required.
 - The usual collection amount is a very small and because it is collected while in lumbar anesthesia so exception care is required.
 - The collected sample should be tested as soon as possible to avoid cell transformation.
1. Right after or within 2hrs of storing refrigerator, put collected sample into 50ml centrifuge tube and centrifuge in 1800rpm for 5 minutes then pour into Cellprep vial for 30mins for fixation.
 - In case of a small amount of collected specimens (1~5ml), directly pour into Cellprep vial as it's not necessary to remove supernatant.
 - * The characteristic is the same with Urine (low mucus or protein content), so it is fine to use same Cellprep vial.
 2. Label the slide and prepare with Cellprep LBC processor.
 - Before insert Cellprep vial to Cellprep LBC processor, vortex for 5~10 seconds or shake 4~5 times

3. Fix the slide directly in 95% ethanol (more than 30mins), then stain with Papanicolaou staining for examination.
- For CSF specimen, it is recommended to use after 12 hours in order to obtain the best quality of cell morphology

8.3. Respiratory tract (sputum, bronchial washing, bronchial brushing)

- Specimen Collection and precautions: Should not be stored in a Refrigerator for more than 2 hours (Before putting Respiratory tract into Cellprep Vial.

□ Sputum:

- It is fine to put sample directly to Cellprep vial,

But, if you take sputum in large cup at first and then if you put collected sample with the part of blood, discoloration and detail part by a sharp scissors and forceps into Cellprep Vial, it gives high quality of result.

- The proper amount of specimen in vial is about 1ml or less.

□ Bronchial washing:

Wash it with saline and put sucked specimen directly into Cellprep Vial.

- In case of bronchial washing specimens, it is recommended to use after 12 hours in order to obtain the best quality of cell morphology.

□ Bronchial brushing:

Brush can be taken out through pipes and then put it directly into Cellprep Vial.

- Pour all collected specimen from brush into Cellprep Vial.

1. Centrifuge the collected specimen at 1800rpm for 5 mins.

- In case of a lot of blood and mucus at sediment, use lysis solution then process next step.

2. Remove supernatant then pour around 1ml of sediment into Cellprep Vial directly and pre-prepare for more than 30 minutes.

3. Label the slide and prepare with Cellprep LBC Processor.

- Before insert Cellprep Vial to Cellprep LBC System, vortex for 10~20 seconds or shake 5~10 times.

4. Fix the slide directly in 95% Ethanol for more than 30 minutes then stain with Papanicolaou staining method for examination.




8.4. Body fluid (Pleural, Ascites, Peritoneal, Pericardial) & Fine Needle Aspiration cytology (Thyroid, Breast, Salivary gland, Lymph node, Lung, Abdominal organs)

Body fluid & Fine Needle Aspiration cytology

- Specimen Collection and Precautions: Should not stored in a Refrigerator more than 12 hours (Before putting the samples into Cellprep Vial.
- Collection: Insert the fine needle to the part where user wants and check the part and fix. Connect Fluid set and drain it to empty bottle.

1. Put the collected specimen into 50ml centrifuge tube then centrifuge in 1800rpm for 5 minutes.
 - But, if the specimen collected is a very small dose (1~2ml) insert it straight into the Cellprep vial without any pre-preparation.
 - In case of a lot of blood and mucus at sediment, use lysis Solution then process next step.
2. Remove supernatant then pour around 1ml of sediment into Cellprep Vial directly and pre-prepare for more than 30 minutes.
3. Label the slide and prepare with Cellprep LBC Processor.
 - Before insert Cellprep Vial to Cellprep LBC System, vortex for 10~20 seconds or shake 5~10 times.
4. Fix the slide directly in 95% Ethanol more than 30 minutes then stain with Papanicolaou staining method for examination.
 - In case of Body Fluid Specimen, a better quality slide is produced if the smear is made after being left to fixate for 12 hours in the Cellprep Vial.

8.5. Additional information for preparing liquid specimens (Body Fluid, CFS, Bronchial washing)

- All the Cellprep Vials can be used after 30minuts from the putting the specimen into the vial.
- If degeneration is shown in the first smear, repeat the specimen after 12 hours.
- Only in case of bronchial washing, CSF, Body fluid, it is recommended to use after 12 hours. It is because of Cellprep's unique components (elements) to prevent fatal degeneration (Shrink, Hyper chromaticness)

Part. 9. Recommended Staining Method of Cellprep

9.1. Recommended Staining Method of Cellprep (manual/Gyn)

Step No.	Name of Chemical for Stain	Duration	
		Gyn 30 minutes ¹	Non Gyn 30 minutes ²
1	80% Ethanol	10 dips	10 dips
2	70% Ethanol	10 dips	10 dips
3	Tap Water	2 minute	2 minute
4	✓ Hematoxylin	4 minute	2 minute
5	Tap Water	1 minute 30 seconds	1 minute 30 seconds
6	0.5% HCL	3~5dips	3~5dips
7	Tap Water	5 minute	5 minute
8	95% Ethanol	10 dips	10 dips
9	95% Ethanol	10 dips	10 dips
10	✓ OG-6	3 minute	1 minute 30 seconds
11	95% Ethanol	10 dips	10 dips
12	95% Ethanol	10 dips	10 dips
13	✓ EA-50	3 minute	1 minute 30 seconds
14	95% Ethanol	10 dips	10 dips
15	95% Ethanol	10 dips	10 dips
16	100% Ethanol	10 dips	10 dips
17	100% Ethanol	10 dips	10 dips
18	Xylene	10 dips	10 dips
19	Xylene	1 minute(10 dips)	1 minute(10 dips)

What reagent user want	Expected manufacture concentration	95% ethanol (ml)	DW (ml)
50% alcohol	50%	526	474
70% alcohol	70%	737	263
80% alcohol	80%	840	160
What reagent user want	Expected manufacture concentration	70% ethanol(ml)	Hcl (ml)
0.5% hcl alcohol	0.5%	995	5

9.2. Recommended Staining Method of Cellprep (Auto Stainer)

1. Fixation step: the slides which inserted in rack have to be fixed for 30 minute in 95% Alcohol.

* Please exchange 'the fixation reagent' for new one everyday

- Wash step by flowing water: wash the fixed specimen with flowing water for 5~7minutes.
- Pap stain (Auto stainer : 1minute each)

[1] harris hematoxylin ①	▶	[2] harris hematoxylin ②	▶	[3] harris hematoxylin ③	▶
[4] Wash step by flowing water	▶	[5] Wash step by flowing water	▶	[6] 0.5% HCL	▶
[7] Wash step by flowing water	▶	[8] Wash step by flowing water	▶	[9] Wash step by flowing water	▶
[10] Wash step by flowing water	▶	[11] 95% Alcohol	▶	[12] 95% Alcohol	▶
[13] 95% Alcohol	▶	[14] OG-6	▶	[15] 95% Alcohol	▶
[16] 95% Alcohol	▶	[17] EA-50 ①	▶	[18] EA-50 ②	▶
[19] 95% Alcohol	▶	[20] 95% Alcohol	▶	[21] 95% Alcohol	▶
[22] 99% Alcohol	▶	[23] 99% Alcohol	▶	[24] Xylene	▶
[25] Xylene					

Pressure & Time


Specimen	Mode	Suction			Blower		
		Pressure	Time(sec)		Pressure	Time(sec)	
GYNE	A	10.0	35.0	1.0	0.05	0.05	0.5
	B	15.0	35.0	1.0	0.05	0.05	0.5
	C	1.3	30.0	1.0	0.05	0.05	0.5
URINE	A	20.0	40.0	1.0	0.05	0.05	0.5
	B	27.0	40.0	1.0	0.05	0.05	0.5
	C	15.0	35.0	1.0	0.05	0.05	0.5
BODY FLUID	A	10.0	35.0	1.0	0.05	0.05	0.5
	B	20.0	40.0	1.0	0.05	0.05	0.5
	C	4.0	35.0	1.0	0.05	0.05	0.5
FNA	A	20.0	40.0	1.0	0.05	0.05	0.5
	B	37.0	40.0	1.0	0.05	0.05	0.5
	C	10.0	35.0	1.0	0.05	0.05	0.5
RPT	A	20.0	35.0	1.0	0.05	0.05	0.5
	B	37.0	35.0	1.0	0.05	0.05	0.5
	C	10.0	35.0	1.0	0.05	0.05	0.5

9. Used symbols and the marking are presented in new Section: "Symbol identification and marking" as follows

Part. 11. Symbol identification and marking

11.1. Symbols identification

Symbol	Name of the symbol	Description
	Caution! 	Indicates that the user should to study the important safety-related information in the User Manual, in particular, warnings and precautions, which couldn't be applied to the medical device itself for various reasons.
	Warning. Danger.	Used to attract the attention to the types of danger. This symbol is used with an additional safety symbol and an explanatory inscription.
	Manufacturer	Indicates the manufacturer of the medical device.
	Declaration of conformity.	CE logo: this instrument complies to European Community directives.
	In-vitro diagnostic.	IVD medical device according to 98/79/EC directive

	Symbol WEEE	WEEE European directive symbol (2002/96/EC): electric/electronic device do not throw out in the environment
	Serial number	The manufacturer can identify a particular medical product by the serial number
	Follow the User Manual	Need to follow recommendation of User Manual
	Highly flammable	This symbol is indicated flammable reagents, solvents and cleaning agents.
	Electrical Hazard	Placed on electrical equipment and instruments, doors of power plates, on electrical panels and cabinets, as well as on fences of current-carrying parts of equipment, mechanisms, instruments.
	ON	By clicking on this sign, the unit is turned on.
	OFF	By clicking on this sign, the unit is turned off.
	Alternative current	Indicates alternative current
	Nominal voltage, V	Indicates nominal voltage
	Max. nominal power, kW	Indicates max. nominal power
	Warning. Moving parts. Hand injury is possible.	This symbol is placed on equipment, equipment assemblies, covers and doors, where it is possible to injure hands.
	UP	Indicates the correct vertical position of equipment.
	Protect from moisture	Necessary to protect the cargo from moisture
	Don't take by hooks	Prohibition of using hooks for lifting instrument
	Fragile, handle with care	Indicates that a medical device could be break or damaged by inappropriate handling.
	Fragile, side impacts are prohibited	Indicates that a medical device could be break or damaged by side impact, if not properly handled.
	Handle with care	The symbol is indicated that the device in the package requires careful handling.

11.2. Marking

NUMBER OF REGISTERED MANUFACTURE
REGISTRATION NUMBER OF PRODUCT 1828
GYEONGIN JESIN 07-39

MODEL NAME CELLPREP AUTO LBC SYSTEM
PRODUCTION NUMBER XX-XXXXXXXX-XXX*
DATE OF PRODUCTION YYYYYYYY*
WEIGHT 175 kg
DIMENSION 1,000mm x 650mm x 800mm
(W) X (D) X (H)
ELECTRIC POWER SOURCE 220-230 V ~, 50/60Hz
POWER CONSUMPTION Max 1 KW
PACKING UNIT 1 DEVICE / SET

USING METHOD AND ATTENTIONS
FUNCTION REFER TO USERS MANUAL
GRADE AND TYPE OF PROTECTION ON ELECTRICSHOCKS
GRADE 1 DEVICE, TYPE B EQUIPMENT

MANUFACTURER:  BIODYNE Co., Ltd.
ADDRESS #301, KRANZ TECHNO Bld., Dunchon-daero 388,
Jungwon-gu, Seongnam-si, Gyeonggi-do,
Korea
TEL / FAX (031) 741-7253 / (031) 741-7259
THIS DEVICE IS MEDICAL EQUIPMENT

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* XX-XXXXXXXX-XXX

YYYYYYYY - variable data

11.3. Warning notices and symbols on the instrument

CAUTION!
HAZARDOUS VOLTAGE.
Contact may cause electric shock or burn.
TURN OFF and LOCK OUT power before servicing.

WARNING!
MOVING PARTS CAN CRUSH AND CUT.
Do not operate with guard removed.
Follow lock out procedure before servicing.

WARNING. DANGER.

10. The Disposal information is presented in new Section: "Disposal information"

Part. 12. Disposal information

12.1. Disposal information of Cellprep AUTO

Before disposing is required to clean and disinfect the instrument.

The Cellprep AUTO and/or its parts and accessories must be disposed in accordance with the existing applicable, local regulations. For the Russian Federation it is recommended to be disposed in accordance with the rules implemented by "СанПиН 2.1.7.2790-10 (Sanitary and epidemiological requirements for the medical waste handling)" for the medical waste referred to class A (non-hazardous waste), excluding electronic parts.

Electronic parts must be disposed in accordance with local specified regulations, in relation to household and industrial waste (in particular, office and computer equipment).

This symbol indicates separate collection
for electrical and electronic equipment.

12.2. Disposal of used materials

Waste solutions, filters and glasses should be treated as biohazardous materials and disposed of in accordance with federal, municipal and local government regulations.

For the Russian Federation it is recommended to be disposed in accordance with the rules implemented by "СанПиН 2.1.7.2790-10 (Sanitary and epidemiological requirements for the medical waste handling)" for the medical waste referred to class Б (epidemiologically hazardous waste).

11. The information about the applicable national standards used by the manufacturer for the medical product is contained in the new section: "List of applicable standards"

13. List of applicable standards

Cellprep AUTO System meets the following standards:

- 98/79/EEC In vitro diagnostic medical devices Directive;
- 2006/95/EC Low Voltage Directive;
- 2004/108/EC Electromagnetic compatibility EMC Directive;
- EN 61326-1:2013 Electrical equipment for measurement, control and laboratory use. EMC requirements. General requirements;
- EN 55011:2009/A1:2010 Industrial, scientific and medical equipment. Radio-frequency disturbance characteristics. Limits and methods of measurement;
- EN 61000-3-2:2014 Electromagnetic compatibility (EMC). Limits. Limits for harmonic current emissions (equipment input current 16 A per phase);
- EN 61000-3-3:2013 Electromagnetic compatibility (EMC). Limits. Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current 16 A per phase and not subject to conditional connection;
- EN 61010-1:2010 Safety requirements for electrical equipment for measurement, control, and laboratory use. General requirements;
- EN 61010-2-081:2002/A1:2003 Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-081: particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes;
- EN 61010-1-101:2002 Safety requirements for electrical equipment for measurement, control and laboratory use. Particular requirements for in vitro diagnostic (IVD) medical equipment.
- IEC 62471:2006 Lamps and lamp systems. Photo biological safety.
- EN IEC 62304:2013 Medical devices. Software. Life cycle processes

The User Manual with this Appendix № 1 is valid for Russian Federation only and will be supplied with **Cellprep AUTO LBC System** to the Russian market.

Approved

(sign, stamp)

CCO

Jamin Lee

20.6.20

Blodyne

THE APPENDIX SHOULD BE SEWED TOGETHER with Operator Manual, THIS PACK SHOULD BE SIGNED/STAMPED and NOTARIZED

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Registered No. 2020 - 6650

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NOTARIAL CERTIFICATE

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As a result of checking at my
office, I have found that the

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attached

Cellprep AUTO LBC System user
manual

COPY

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exactly corresponds with the
original.

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This is hereby attested on
this 3rd day of June 2020

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Seoul Central District Prosecutors' Service
(Diplomatic Center B/D 2F, Seocho-dong) 202,
2558, Nambusunhwan-ro, Seocho-gu, Seoul, Korea
Signature of the Notary Public

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C. G. Kim.

KIM CHEOL GEE

This office has been authorized by
the Minister of Justice, the Republic
of Korea to act as Notary Public
since March 14, 2007
under Law No. 7428.

APOSTILLE

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1. Country : Republic of Korea

This public document

2. has been signed by KIM CHEOL GEE

3. acting in the capacity of Notary Public

4. bears the seal/stamp of HANMI LAW AND NOTARY OFFICE

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