65533

JD

APPROVED

Reliance Orthodontic Products, Inc.

Organization

1540 WEST THORNDALE AVE.
ITASCA, IL USA 60143

Address

QA/RA Specialist

Position

Sylvia Trenkler

Name

02.11.2021

Date DD.MM.YYYY

Sylvia Trenkler

Signature

**Instructions for Use
for medical device:
Orthodontic Products Bonding Materials in different variants, with
accessories**

State of Illinois

County of Cook

Signed and attested before me on

Nov. 2, 2021 (Date)

By Sylvia Trenkler (Name of Person)

Signature of Notary Public

2021

Instructions for use
Orthodontic Products Bonding Materials in different variants, with accessories
Orthodontic Products Bonding Materials

Perfect-A-Smile Pontic Aligner Paint

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Orthodontic Products Bonding Materials in different variants, with accessories

I. Orthodontic Products Bonding Materials

1. Perfect-A-Smile; PP/Bleach shade (catalogue number PPBL):
 - Perfect-A-Smile; PP/Bleach shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
2. Perfect-A-Smile; PP/A1 shade (catalogue number PPA1):
 - Perfect-A-Smile; PP/A1 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
3. Perfect-A-Smile; PP/A2 shade (catalogue number PPA2):
 - Perfect-A-Smile; PP/A2 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
4. Perfect-A-Smile; PP/A3.5 shade (catalogue number PPA35):
 - Perfect-A-Smile; PP/ A3.5 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
5. Perfect-A-Smile; PP/B1 shade (catalogue number PPB1):
 - Perfect-A-Smile; PP/B1 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
6. Perfect-A-Smile; PP/C1 shade (catalogue number PPC1):
 - Perfect-A-Smile; PP/C1 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
7. Perfect-A-Smile; PP/C4 shade (catalogue number PPC4):
 - Perfect-A-Smile; PP/C4 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
8. Perfect-A-Smile; PP/D3 shade (catalogue number PPD3):
 - Perfect-A-Smile; PP/D3 shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
9. Perfect-A-Smile; PP/Light Shade (catalogue number PASL):
 - Perfect-A-Smile; PP/Light Shade, 4g jar – 1 pc.;
 - spatula - 1 pc.;
 - micro-brushes - 10 pcs. (if needed);
 - instructions for use – 1 pc.
10. Perfect-A-Smile; PP/Medium Shade (catalogue number PASM):
 - Perfect-A-Smile; PP/Medium Shade, 4g jar – 1 pc.;

- spatula - 1 pc.;
- micro-brushes - 10 pcs. (if needed);
- instructions for use – 1 pc.

Further in the text can appear as - "product", " medical device"

Product Description

The Perfect·A·Smile Pontic Paint is intended for use as a light cure coloring agent of thermoplastic aligners.

Potential Users

Trained dental professionals.

Indications, Contraindications and Adverse events

- Indications

- bracket adhesion.
- bonding lingual retainers and splinting materials,
- creation of archwire stops
- indirect bonding of orthodontic brackets and appliances.
- retention of thermoplastic aligners
- creation of occlusal buildups

- Contraindications

There are no known contraindications.

- Adverse events (Potential complications)

Skin sensation in susceptible persons (Individual intolerance to product components).


Indications for use: The Perfect·A·Smile™ Pontic Paint is intended for use as a light cure coloring agent of thermoplastic aligners.

This patented formula works with clear thermoforming aligner plastics.

Contraindication: Perfect A Smile is not recommended for use with cloudy, (C) thermoforming materials.

Instructions for use:

1. To ensure uniformity due to settlement of the product, mix the pontic paint thoroughly with a spatula until uniform. **PRODUCT MUST BE MIXED BEFORE EACH USE.**
2. Remove the amount desired and place on mixing pad. **IMMEDIATELY REPLACE THE LID TO PROTECT FROM LIGHT.**

3. Apply 3-4 coats of the pontic paint with a microbrush applicator, stroking over each coat several times. Light cure for 10 seconds at close range with a dental curing light after each coat. **NOTE: Shades can be mixed ON A MIXING PAD to achieve the color desired.**

4. After final coat is applied and light cured, perform a final light cure through the aligner for 20 seconds from the labial and lingual surfaces.

Please Note: For the most accurate matching with your aligner material, fabricate a shade guide by placing 3 coats of each shade on a wax covered mixing pad or an any clear aligner plastic. Light cure for 10 seconds between coats. If using a mixing pad, the final polymerized template can then be lifted off the surface and held up to the tooth for matching purposes.


Precautions: Unpolymerized resins may cause skin sensation in susceptible persons. In case of contact with the skin, wash thoroughly with soap and water. Contact with other tissues: rinse spill area immediately with copious amounts of water for several minutes. **DO NOT INHALE FUMES. USE IN A WELL VENTILATED AREA.**

MicroBrush applicators are single use devices. It is required that a new microbrush applicator is used for each application.

Use of UV, halogen and LED curing lights can be associated with eye damage. The patient, doctor and the assistant should wear properly designed protective eyewear when using these light curing devices.

Creating a temporary tooth on a stone model

Replace the missing tooth or teeth with any denture tooth shaped to properly fill the edentulous space. Flowtain™ can be used to hold the denture tooth in place.

Roll approximately ½ gram of Light Bond™ Operative between your index and thumb.

Place roll of composite on gingival surface of the missing tooth socket and lightly depress ball of composite between thumb and index to produce incisal height slightly shorter than adjacent teeth*. Precise shaping of the incisal, mesial and distal can be done with an instrument.

Light Cure for 30 seconds from 2 (two) angles

Fractured or worn teeth can be repaired or altered in shape with Light Bond Operative.

***IMPORTANT!** For a single pontic tooth on the model, create the incisal edge slightly shorter (approx. 0.5mm) than the adjacent teeth. This compensates for the thickness of the aligner material on the facial surface of the adjacent teeth once the aligner socket is painted.

For multiple or adjacent teeth, (for an exact match), form the pontic tooth on the model to the exact incisal height of adjacent teeth. To allow room for the paint on existing teeth in the aligner socket, apply a thin coat of FlowTain™ on the facial surfaces of adjacent teeth - similar to the protocol for bleaching relief voids in bleaching trays.

Storage and Stability

Store at room temperature (15-30 °C) and humidity under 80%. Allow product to reach room temperature prior to use.

List of pharmaceutical substances

Not applicable. Products do not contain any pharmaceutical substances

List of Materials of Animal and (or) Human Origin

Not applicable. Products do not contain any materials of human/animal origin

Disposal

Disposal methods: Dispose of content/container in accordance with the local/regional/national/international/regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organisation and implementation of sanitary and anti-epidemic (preventive) measures" to the waste class B.

Description of the sterilization method, information about validation methods for the sterilization process, and validation of the packaging process for sterile products

Sterilization is not applicable. These products are not supplied sterile and do not require sterilization.

Warranty

Reliance Orthodontic Products, Inc. recognizes its responsibility to replace products if proven to be defective. Reliance Orthodontic Products, Inc. does not accept liability for any damages or loss, direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connections therewith.

Package content

For Perfect-A-Smile; PP/Bleach shade (catalogue number PPBL):

- Perfect-A-Smile; PP/Bleach shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/A1 shade (catalogue number PPA1):

- Perfect-A-Smile; PP/A1 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/A2 shade (catalogue number PPA2):

- Perfect-A-Smile; PP/A2 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/A3.5 shade (catalogue number PPA35):

- Perfect-A-Smile; PP/ A3.5 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/B1 shade (catalogue number PPB1):

- Perfect-A-Smile; PP/B1 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/C1 shade (catalogue number PPC1):

- Perfect-A-Smile; PP/C1 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/C4 shade (catalogue number PPC4):

- Perfect-A-Smile; PP/C4 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/D3 shade (catalogue number PPD3):

- Perfect-A-Smile; PP/D3 shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/Light Shade (catalogue number PASL):

- Perfect-A-Smile; PP/Light Shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

For Perfect-A-Smile; PP/Medium Shade (catalogue number PASM):

- Perfect-A-Smile; PP/Medium Shade, 4g jar – 1 pc.;
- spatula - 1 pc.;
- micro-brushes - 10 pcs.;
- instructions for use – 1 pc.

Manufacturer of the medical device

Developer

Reliance Orthodontic Products, Inc.
1540 WEST THORNDALE AVE. Itasca, IL USA 60143

Manufacturer

Reliance Orthodontic Products, Inc.
1540 WEST THORNDALE AVE. Itasca, IL USA 60143

Addresses of Manufacturing Sites for medical device

Reliance Orthodontic Products, Inc.
1540 WEST THORNDALE AVE. Itasca, IL USA 60143

Customer Information

Rx Only: Federal law restricts this device to sale or on the order of a dental professional.


Customer Service

Authorized representative of the manufacturer in the Russian Federation:

“Ortholight MSK” LLC, 5 Golovinskoe highway, Building 1, Floor 9, room 9024, Moscow, 125212, Russia. Phone +7 (495) 136-76-76

Symbols used on labels and packaging

	Use by date		Do not reuse
	Catalogue number		Lot number
	Caution	CONTENTS	Contained items

	Conformity sign		Product meets the requirements of the applicable EC directives
	European Authorized Representative		
	Manufacturer		Temperature limit (storage conditions)

	Recycling code (Non-corrugated fiberboard)		Recycling code (Low-density polyethylene)
	Green dot		

	Flammable		Corrosive substance
--	-----------	---	---------------------

	Poison (Skull and crossbones)		Corrosive substance
	Harmful		

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